



1076214

INSTRUCTIONS

This spreadsheet is designed to record the raw fiber counts for air and dust samples analyzed by TEM

This is version TEM32d


Raw Data Recording

Raw data are to be recorded by the analyst in hard copy using Lab sheet 1 and as many Lab Continuation pages as needed

Data Entry

Data on the hard copy lab sheets are to be entered electronically on data entry sheets 1 and 2

Areas for data entry are highlighted in **YELLOW**



OR are indicated by a **PULL-DOWN MENU**



Cells that are shaded gray do not require any data input

Cells that are shaded red either require data input or contain an apparent data inconsistency

Do not enter data in any other location!

Avoid the drag and drop method for copying.

Enter all values individually

File Save

After entering all data on Data Entry 1 (or any time thereafter), SAVE THE FILE by clicking on the macro button located on data Entry 1 or Data Entry 2.

The file name is generated automatically by concatenating information provided in Data Entry 1, in

a) EPA Index ID, b) Lab Sample Number, c) Counting Rule, d) Analysis Date, e) Prep method, and for example, if the EPA number is **2-012345** and the lab ID number is **abc-1234**, the method is

AHERA, the analysis date is **03-04-01**, the prep method is **Direct**, and the QA type is **Recount Same**

the file name will be: "**2-012345_abc-1234_AHERA_03-04-01_DRS.XLS**"

opened:

--if you open EXCEL and then open the TEM EDD, the new file will be saved in the same directory as the TEM EDD.

--if you open the TEM EDD from Explorer, then the file will be saved in you default directory for Excel (usually this is C:\Documents and Settings\My Documents)

TEM - General Counting Rules

Follow established project-specific counting and stopping rules.
Do not count non-asbestos material (NAM) structures.

ISO 10312 Counting Rules:

The lab analyst is responsible for identifying if a Bundle (B) is countable (i.e., substructures within have an aspect ratio greater than or equal to the specified aspect ratio). Aspect ratio rules do not apply to compact clusters, compact matrices, or residuals.
If a structure is non-countable (e.g., crosses the top or left grid bar), it should be identified with a "0" in the Total column.
All countable structures should be identified with a non-zero sequential number in the Total column. Enter a non-zero sequential number for any unique primary structure in the first row only. If there are multiple rows associated with the primary structure, it is not necessary to repeat the primary number in every row.
NOTE: Data Entry 2, column R is designed to flag any potential data entry errors associated with countable/non-countable assignments.

AHERA/ASTM Counting Rules:

Valid Structure Types are Fiber (F), Matrix (M), Cluster (C), and Bundle (B).

For Matrices (M), record the dimensions of the protrusion, not the matrix.

Note: Currently, the spreadsheet will incorrectly exclude one particular class of matrix structure from the AHERA binned results summary -- a matrix for which a protruding bundle has an aspect ratio less than 5:1, but the underlying fibrils meet the aspect ratio requirement. This is a limitation of the TEM spreadsheet.

If this type of structure is encountered, the analyst should note this issue in the analysis comments.

The lab analyst is responsible for identifying if a Cluster (C) or Bundle (B) is countable (i.e., substructures within have an aspect ratio greater than the specified aspect ratio). Do not record non-countable structures.

All countable structures should be identified with a non-zero sequential number in the Total column.

NOTE: Data Entry 3, column R is designed to flag any potential data entry errors associated with countable/non-countable assignments.

DATA ENTRY INSTRUCTIONS

| Data Item | Description | Notes |
|---|---|---|
| DATA ENTRY 1 | | |
| File Type | Correction/revision status of the file | Select the appropriate File Type (original, correction 1, correction 2, correction 3, correction 4) from the pull-down list. The filename will be automatically updated to reflect the correction status (i.e., ".C1" suffix will be added to the "correction 1" file). It is not necessary to manually change the filename. |
| Laboratory Name | Name of lab performing analysis | Use a standard name for all sheets. |
| Instrument | Instrument used for analysis | |
| Voltage (KV) | Voltage used for analysis | |
| Magnification | Magnification used for analysis | |
| Grid opening area | Size of grid opening | Enter as a value in units of mm ² . This field should not contain any text. |
| Scale: 1L = | Scale adjustment factor for length | Enter 1.0 if dimensions are expressed as um. Otherwise, enter um per unit screen length. |
| Scale: 1D = | Scale adjustment factor for width | Enter 1.0 if dimensions are expressed as um. Otherwise, enter um per unit screen width. |
| Primary filter area | Size of primary effective filter area (EFA) | Enter as a value in units of mm ² . This field should not contain any text. |
| Secondary filter area | Size of secondary effective filter area (EFA) | Enter as a value in units of mm ² . This field should not contain any text. |
| Category | Sample type | Select the appropriate category (Field, Replicate, Duplicate, Blank) from the pull-down list. |
| Filter Status | Status of the sample filter | Select the appropriate category (Analyzed, Overloaded, Damaged, Missing, Cancelled) from the pull-down list. |
| EPA Index ID | Unique sample identifier | Enter the EPA Index ID exactly as it appears on the sample (ie: 2-00013). This field should not contain any spaces or a lab QA type suffix. |
| Sample Type | Sample media | Select the appropriate category (Air, Dust, Dustfall) from the pull-down list. |
| Air volume (L), dust area (cm2), or dustfall container area (cm2) | Air volume (L), dust area (cm2), or dustfall container area (cm2) | If the sample media is air, enter the air volume in units of L. If the sample media is dust, enter the dust sample area in units of cm ² . If the sample media is dustfall, enter the dustfall container area in units of cm ² . This field should not contain any text. |
| Date Received | Date sample was received by lab | Enter as a valid date with the format MM/DD/YYYY. |
| Lab Job Number | Job number assigned to analysis by lab | |
| Lab Sample Number | Sample identifier as designated by lab | If the sample is a lab QA Type, DO NOT add the QA suffix to the Lab Sample ID. If the file has been corrected or modified, DO NOT add "REV" to the Lab Sample ID. Use the File Type pull-down list to note the revision status. |
| Number of grids prepared | Number of grids prepared by lab | |
| Prepared by | Name of lab preparation personnel | Enter as first initial and last name (J. Smith) |
| Preparation date | Date sample was prepared | Enter as a valid date with the format MM/DD/YYYY. |
| EPA COC Number | EPA chain of custody number | |
| Analyzed by | Name of analyst | Enter as first initial and last name (J. Smith) |
| Date Analyzed | Date sample was analyzed by lab | Enter as a valid date with the format MM/DD/YYYY. |
| Prep | Analysis Prep Type | Select the appropriate Prep type (Direct, Indirect, Indirect - Ashed) from the pull-down list. |
| Loose Material in the Cow? | Is there loose material or debris in the saw? | If the sample media is air, select either "yes" or "no" from the pull-down list. |
| Counting Rules | Counting Rules utilized to analyze sample | Select the appropriate Counting Rule from the pull-down list. If the sample media is dust, select either ASTM-AHERA or ISO. If the sample media is air, select either AHERA or ISO. |
| Grid Storage Location | Grid storage location identifier | |
| F-factor Calculation: | | For Direct prep, F-factor is set to equal 1 and specific F-factor calculations are not required. |
| Indirect Prep inputs: | | |
| Fraction of primary filter used | | Enter the fraction as a value. This field should not contain any text. |
| First resuspension volume or rinaste volume (mL) | | Enter as a value in units of mL. This field should not contain any text. |
| Volume applied to secondary filter (mL) or used for serial dilution | | Enter as a value in units of mL. This field should not contain any text. |
| Inputs for Serial Dilutions: | | |
| Second resuspension volume (mL) | | Enter as a value in units of mL. This field should not contain any text. |
| Volume applied to secondary filter (mL) or used for serial dilution | | Enter as a value in units of mL. This field should not contain any text. |
| Third resuspension volume (mL) | | Enter as a value in units of mL. This field should not contain any text. |
| Volume applied to secondary filter (mL) or used for serial dilution | | Enter as a value in units of mL. This field should not contain any text. |
| Inputs for Ashing of Secondary Filter: (Note: only complete this cell if the SECONDARY filter was ashed. Leave blank if the primary filter was ashed.) | | |
| Fraction of secondary filter used | | Enter the fraction as a value. This field should not contain any text. |
| QA Type | Lab quality control code | Select the appropriate laboratory QA type (Not QA, Recount Same, Recount Different, Repreparation, Verified Analysis, Reconciliation, Lab Blank, Interlab) from the pull-down list. |
| Comments | Sample/Analysis comments | |
| Recording Rules: | | |
| Minimum Aspect Ratio | Minimum aspect ratio of structures that are to be recorded | Select the appropriate minimum aspect ratio (none, ≥ 3:1, ≥ 5:1) from the pull-down list. |
| Minimum Length | Minimum length of structures that are to be recorded | Enter as a value in units of um. If no minimum length is specified, enter 0.00. |
| Minimum Width | Minimum width of structures that are to be recorded | Enter as a value in units of um. If no minimum length is specified, enter 0.00. |
| Stopping Rules (stop if any one of the following criteria are met): | | |
| Target Sensitivity | Target sensitivity | Enter as a value in units of slices for air or skin* for dust or dustfall. This field should not contain any text. If no target sensitivity is required, leave this field blank. |
| Maximum # of GOs | Maximum # of GOs | Enter as a value. This field should not contain any text. If no maximum # of GOs has been established, leave this field blank. |
| Maximum # of Structures | Target # of structures | Enter as a value. This field should not contain any text. If no target # of structures is required, leave this field blank. |
| DATA ENTRY 2 | | |
| Data Entry by | Name of data entry personnel | Enter as first initial and last name (J. Smith) |
| Data Entry Date | Date results were entered | Enter as a valid date with the format MM/DD/YYYY. |
| QA by | Name of QA personnel | Enter as first initial and last name (J. Smith) |
| QA Date | Date results were QA'd | Enter as a valid date with the format MM/DD/YYYY. |
| Grid | Grid identifier | Enter the appropriate grid in the first row only. If there are multiple rows associated with the grid, it is not necessary to repeat the grid in every row. |
| Grid Opening | Grid opening location identifier | Enter the appropriate grid opening in the first row only. If there are multiple rows associated with the grid opening, it is not necessary to repeat the grid opening in every row. Indicate grid openings that DO NOT count chrysotile with an "x" following the grid opening name (eg. A-4"). Indicate grid openings that COUNTED ONLY A FRACTION OF THE GO for chrysotile with an "f" following the grid opening name (eg. A-4f). If only a fraction of the total GO was counted for chrysotile, you must enter the fraction of the GO counted in the appropriate column (Fract. GO Chrys.). |
| Structure Type | Structure Type code | Enter the structure type code. For ISO, see the analytical method for valid structure type codes. For AHERA, valid structure types are Fiber (F), Matrix (M), Cluster (C), and Bundle (B). |
| Primary | Primary Structure identifier | For ISO, enter a non-zero sequential number for any unique primary structure in the first row only. If there are multiple rows associated with the primary structure, it is not necessary to repeat the primary number in every row. This field is not used for AHERA or ASTM. |
| Total | Total Structure identifier | Assign a "0" to any non-countable or excluded structure and a non-zero number to any countable structure. **See specific instructions above for details on populating this field for ISO and AHERA/ASTM. |
| Length | Structure length | Enter dimensions either in absolute units (um) or in screen units. |
| Width | Structure width | Enter dimensions either in absolute units (um) or in screen units. |
| Identification | Identification | |
| Mineral Class | Description of the structure mineral class type | Enter a "1" in the appropriate column; choices include "LA" (Libby amphibole), "OA" (other amphibole), "C" (chrysotile). You may only select one mineral class type for each structure. NOTE: Do not enter "NAM" (non-asbestos material) structures. |
| Comments | Structure comments | |
| Sketch | Sketch | Enter a "1" in this column if yes, otherwise leave this field blank. |
| Photo | Photo | Enter a "1" in this column if yes, otherwise leave this field blank. |
| EDS | EDS | Enter a "1" in this column if yes, otherwise leave this field blank. |
| Fract. GO Chrys. | Fraction of GO counted for Chrysotile | If only a fraction of the GO was counted for chrysotile, enter the fraction counted as a value between 0 and 1 (eg: if only 1/4 of the GO was counted for chrysotile, the value entered would be 0.25). |

Enter Project Name Here:

LIBBY

LIBBY
TEM Asbestos Structure Count

| | |
|--|-------------|
| Laboratory name: | RESI |
| Instrument | JEOL 100 CX |
| Voltage (KV) | 100 KV |
| Magnification | 20,000 |
| Grid opening area (mm ²) | 0.011 |
| Scale: 1L = | 0.28 um |
| Scale: 1D = | 0.056 um |
| Primary filter area (mm ²) | 385 |
| Secondary Filter Area (mm ²) | 201 |
| Category (Field, Rep., Dup., Blank) | Field |
| Primary filter pore size (um) | 0.45 |

| | |
|---|------|
| EPA Sample Number: | |
| Sample Type (A=Air, D=Dust, DF = Dustfall): | |
| Air volume (L), dust area (cm ²), or dustfall container area (cm ²) | |
| Date received by lab | |
| Lab Job Number: | |
| Lab Sample Number: | |
| Number of grids prepared | 3 |
| Prepared by | |
| Preparation date | |
| EPA COC Number: | |
| Secondary filter pore size (um) | 0.22 |

| | |
|--|--------|
| Analyzed by | |
| Analysis date | |
| Method (D=Direct, I=Indirect, IA=Indirect, ashed) | |
| If sample type = air, is there loose material or debris in the cowl? (Yes, No) | |
| Counting rules (ISO, AHERA, ASTM) | |
| Grid storage location | |
| Archive filter(s) storage location | |
| QA Type (Not QA, Recount Same, Recount Different, Re-prep, Verified Analysis, Reconciliation, Lab Blank, Interlab) | Not QA |

F-Factor Calculation (Indirect Preps Only):

Enter data in appropriate cells provided to the right----->

| | | |
|------------------------------------|-------|-------|
| <u>Recording Rules:</u> | | |
| Minimum Aspect Ratio (circle one): | | |
| none | ≥ 3:1 | ≥ 5:1 |
| Minimum Length (um): _____ | | |
| Minimum Width (um): _____ | | |

| | |
|-------------------------|-------|
| <u>Stopping Rules:</u> | |
| Target Sensitivity: | _____ |
| Max # of GOs: | _____ |
| Target # of Structures: | _____ |

| Grid | Grid Opening | Structure Type | No. of Structures | | Dimensions | | Identification | Mineral Class (see below) | | | | Sketch/ Comments | 1 = yes, blank = no | | | Fract. GO Chrys. |
|------|--------------|----------------|-------------------|-------|------------|-------|----------------|---------------------------|----|---|-----|------------------|---------------------|-------|-----|------------------|
| | | | Primary | Total | Length | Width | | LA | OA | C | NAM | | Sketch | Photo | EDS | |
| | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | |

F-factor Calculation:

Indirect Prep Inputs

| | |
|--|--|
| | Fraction of primary filter used for indirect prep or ashing <i>[For dust and dustfall, enter 1.0]</i> |
| | First resuspension volume or rinsate volume (mL) |
| | Volume applied to secondary filter (mL) or used for serial dilution |

Inputs for Serial Dilutions

| | |
|--|---|
| | Second resuspension volume (mL) |
| | Volume applied to secondary filter (mL) or used for serial dilution |
| | Third resuspension volume (mL) |
| | Volume applied to secondary filter (mL) |

Input for Ashing of Secondary Filter

| | |
|--|--|
| | Fraction of secondary filter used for ashing |
|--|--|

LA = Libby-type amphibole

OA = Other (non-Libby type) amphibole

C = Chrysotile

NAM = Non-asbestos material

| |
|--|
| |
|--|

| |
|--|
| |
|--|

| |
|--|
| |
|--|

| |
|--|
| |
|--|

Page 10 of 10

[illegible]

STRUCTURE LOCATIONS WITHIN GRID OPENING

Lab Name: _____ Lab Job Number: _____

Index ID: _____ Lab Sample ID: _____

Lab QC Type (circle one): Reprep for interlab Interlab

Grid: _____ Grid Opening: _____

upper
left
corner

traverse direction

A large empty rectangular box for drawing or sketching structure locations within the grid opening.

Comments:

A rectangular box for writing comments.

FILE NAME: FC-00007 143576-176058 ISO 10-20-07 IA.xls

FILE TYPE: Original

LIBBY

TEM Asbestos Structure Count

| | |
|-----------------------------|--------------|
| Laboratory name: | RESI |
| Instrument | JEOL 100 C X |
| Voltage (KV) | 100KV |
| Magnification | 20,000 |
| Grid opening area (mm2) | 0.0110 |
| Scale: 1L = | 0.280 |
| Scale: 1D = | 0.056 |
| Primary filter area (mm2) | 385.0 |
| Secondary Filter Area (mm2) | 346.0 |
| Category | Field ▼ |
| Filter Status | Analyzed ▼ |

| | |
|--|---------------|
| EPA Sample Number: | FC-00007 |
| Sample Type | Air ▼ |
| Air volume (L), dust sample area (cm2), or dustfall container area (cm2) | 1200 |
| Date received by lab | 8/27/2007 |
| Lab Job Number: | 143576 |
| Lab Sample Number: | 143576-176058 |
| Number of grids prepared | 3 |
| Prepared by | ND/KS/CL |
| Preparation date | 9/20/2007 |
| EPA COC Number | L12730 |

| | |
|--|---------------------|
| Analyzed by | A. Heitger |
| Analysis date | 10/20/2007 |
| Prep | Indirect - Ashed ▼ |
| If sample type = air, is there loose material or debris in the cowl? | Yes ▼ |
| Counting rules | ISO (Air or Dust) ▼ |
| Grid storage location | 143576 |
| Archive filter(s) storage location | 143576 |
| Enter the appropriate data in the cells to the right to calculate the F-factor >>>>> | |
| F- factor | 0.1 |
| QA Type | Not QA ▼ |

Recording Rules:

| | |
|---------|----------------------|
| ≥ 3:1 ▼ | Minimum Aspect Ratio |
| 0.50 | Minimum Length (um) |
| 0.00 | Minimum Width (um) |

Stopping Rules:

| | |
|---------|--|
| 0.00100 | Target Sensitivity |
| 263 | GOs required to reach target sensitivity |
| 100 | Maximum # of GOs |
| 50 | Maximum # of Structure |
| 100 | Estimated # of GOs |

F-factor Calculation:

Indirect Prep Inputs

| | |
|-----|---|
| 1 | Fraction of primary filter used for indirect prep or ashing [For dust and dustfall, enter 1.0] |
| 100 | First resuspension volume or rinsate volume (mL) |
| 10 | Volume applied to secondary filter (mL) or used for serial dilution |

Inputs for Serial Dilutions

| | |
|--|---|
| | Second resuspension volume (mL) |
| | Volume applied to secondary filter (mL) or used for serial dilution |
| | Third resuspension volume (mL) |
| | Volume applied to secondary filter (mL) or used for serial dilution |

Input for Ashing of Secondary Filter

| | |
|---|--|
| 1 | Fraction of secondary filter used for ashing |
|---|--|

| | |
|-----|----------|
| 0.1 | F-factor |
|-----|----------|

COMMENTS

LIBBY
TEM Asbestos Structure Count

EPA SAMPLE ID:
LAB SAMPLE ID:

| |
|---------------|
| FC-00007 |
| 143576-176058 |

| | |
|-------------|-------|
| Sample Type | Air |
| Count Rule | 10312 |

| | |
|---------|------------------|
| Prep | Indirect - Ashed |
| QA Type | Not QA |

ERROR CHECK

ERROR - Check Values

| | |
|------------------|-------------|
| Data Entry by: | S. Gardalen |
| Data Entry date: | 11/7/2007 |

| | |
|----------|--------------|
| QA by: | G. Vettraino |
| QA date: | 11/7/2007 |

Maximum # of Structures Reached-Complete current GO, then stop.

[illegible]

| | | | | | | | | | | | | | | | | |
|---|------|------|----|----|------|----|-----|---|--|--|---|-----------------|---|---|---|--|
| | | MF | | 19 | 12 | 4 | ADX | 1 | | | | WRTA, NaK | 1 | | | |
| | | B | 20 | 20 | 4 | 4 | ADX | 1 | | | | AT, XX | 1 | | | |
| | | F | 21 | 21 | 8 | 9 | ADX | 1 | | | | WRTA, NaK | 1 | | 1 | |
| | | F | 22 | 22 | 4 | 5 | ADX | 1 | | | | WRTA, NaK | 1 | | | |
| | | F | 0 | 0 | 2.5 | 5 | ADX | | | | 1 | High Ca | 1 | 1 | 1 | |
| | E5-4 | MD10 | 23 | | | | | | | | | | | | | |
| | | MF | | 23 | 15 | 11 | ADX | 1 | | | | AT, XX | 1 | | | |
| | | MD10 | 24 | | | | | | | | | | | | | |
| | | MF | | 24 | 5 | 5 | ADX | 1 | | | | WRTA, NaK | 1 | | | |
| | | B | 25 | 25 | 25 | 19 | ADX | 1 | | | | WRTA, NaK | 1 | | | |
| | | MD11 | 26 | | | | | | | | | | | | | |
| | | MF | | 26 | 20 | 3 | ADX | 1 | | | | WRTA, NaK | 1 | | | |
| | | F | 27 | 27 | 12.5 | 6 | ADX | 1 | | | | WRTA, NaK | 1 | | | |
| B | G3-3 | B | 28 | 28 | 4.5 | 6 | ADX | 1 | | | | WRTA, NaK | 1 | | | |
| | | B | 29 | 29 | 25 | 18 | ADX | 1 | | | | WRTA, NaK | 1 | | | |
| | F3-3 | MD10 | 30 | | | | | | | | | | | | | |
| | | MB | | 30 | 8 | 5 | ADX | 1 | | | | WRTA, NaK | 1 | | | |
| | C3-6 | MD10 | 0 | | | | | | | | | | | | | |
| | | MF | | 0 | 4 | 7 | ADX | 1 | | | | WRTA, NaK <3:1 | 1 | | | |
| | | F | 31 | 31 | 19 | 7 | ADX | 1 | | | | WRTA, NaK | 1 | 1 | | |
| | G4-1 | B | 32 | 32 | 18 | 29 | ADX | 1 | | | | AT, XX | 1 | | 1 | |
| | | F | 33 | 33 | 9 | 4 | ADX | 1 | | | | WRTA, NaK | 1 | | | |
| | | MD11 | 34 | | | | | | | | | | | | | |
| | | MB | | 34 | 21 | 13 | ADX | 1 | | | | AT, XK | 1 | | | |
| | | MD11 | 35 | | | | | | | | | | | | | |
| | | MF | | 35 | 15.5 | 7 | ADX | 1 | | | | AT, XK | 1 | | | |
| | | F | 36 | 36 | 5.5 | 3 | ADX | 1 | | | | WRTA, NaK | 1 | | | |
| | | F | 37 | 37 | 7 | 8 | ADX | 1 | | | | WRTA, NaK | 1 | | | |
| | F4-1 | F | 38 | 38 | 6 | 10 | ADX | 1 | | | | WRTA, XK | 1 | | | |
| | | MF11 | 39 | | | | | | | | | | | | | |
| | | MF | | 39 | 21 | 8 | ADX | 1 | | | | WRTA, NaK | 1 | | | |
| | | MF10 | 40 | | | | | | | | | | | | | |
| | | MF | | 40 | 6 | 7 | ADX | 1 | | | | WRTA, NaK | 1 | | | |
| | | MD11 | 41 | | | | | | | | | | | | | |
| | | MF | | 41 | 62 | 5 | ADX | 1 | | | | WRTA, NaK | 1 | | | |
| | | F | 42 | 42 | 8 | 2 | ADX | 1 | | | | | 1 | | | |
| | | MD10 | 43 | | | | | | | | | | | | | |
| | | MF | | 43 | 14 | 14 | ADX | 1 | | | | WRTA, NaK | 1 | | | |
| | G5-3 | MD10 | 44 | | | | | | | | | | | | | |
| | | MF | | 44 | 17 | 3 | ADX | 1 | | | | WRTA, NaK | 1 | | | |
| | | F | 45 | 45 | 7 | 1 | ADX | 1 | | | | WRTA, NaK | 1 | | | |
| | | MD10 | 46 | | | | | | | | | | | | | |
| | | MF | | 46 | 13.5 | 2 | ADX | 1 | | | | WRTA, NaK | 1 | | | |
| | | F | 0 | 0 | 5.5 | 13 | ADX | 1 | | | | WRTA, NaK < 3:1 | 1 | | | |
| | | B | 47 | 47 | 14 | 15 | ADX | 1 | | | | WRTA, NaK | 1 | | | |
| | | MD11 | 48 | | | | | | | | | | | | | |
| | | MF | | 48 | 61 | 5 | ADX | 1 | | | | WRTA, NaK | 1 | | | |

[illegible]

[illegible]

LIBBY**TEM Asbestos Structure Count -- Not AHERA****SAMPLE ID**

| | |
|-------------------|------------------|
| Status | Analyzed |
| EPA Sample Number | FC-00007 |
| QA Sample Type | Not QA |
| Lab Sample Number | 143576-176058 |
| Sample Type | Air |
| Category | Field |
| Prep | Indirect - Ashed |
| Counting rules | 10312 |

PARAMETERS

| | |
|---------------------------------------|----------|
| Number of Grid Openings (amphibole) | 12 |
| Number of Grid Openings (chrysotile) | 12 |
| Area of grid opening (mm2) | 0.011 |
| F Factor | 1.00E-01 |
| Effective primary filter area (mm2) | 385.0 |
| Effective secondary filter area (mm2) | 346.0 |
| Volume (L) or Sample Area (cm2) | 1200 |
| Area counted (mm2) for LA/OA | 0.132 |
| Area counted (mm2) for C | 0.132 |

Recording Rules:

| | | |
|---------|-----------------|----------------|
| Min. AR | Min length (um) | Min width (um) |
| ≥ 3:1 | 0.5 | 0 |

Stopping Rules:

| | | |
|----------|---------|-------|
| Target S | Max GOs | Max N |
| 0.001 | 263 | 50 |

TOTAL COUNTS

| | |
|---|-----------|
| Total Structures | Not AHERA |
| Total Asbestos structures | Not AHERA |
| Total NAM | Not AHERA |
| Total Not AHERA structures | Not AHERA |
| Total Non-Not AHERA asbestos structures | Not AHERA |

COUNTS (based on countable Not AHERA structures only)

| | | | | |
|-------------------------------|----|----|---|-------|
| | LA | OA | C | Total |
| Not AHERA Structures (< 5 um) | 32 | 0 | 2 | 34 |
| Not AHERA Structures (≥ 5 um) | 19 | 0 | 0 | 19 |

AIR CONCENTRATION (s/cc)

| | | | | | |
|-----------------------------------|---------|----------|---------|---------|------------------------------------|
| | LA | OA | C | Total | |
| Loading on primary filter (s/mm2) | 3.5E+03 | <6.8E+01 | 1.4E+02 | 3.6E+03 | (if Total = ND, DL is based on LA) |
| Air Conc (s/cc) | 1.1E+00 | <2.2E-02 | 4.4E-02 | 1.2E+00 | (if Total = ND, DL is based on LA) |

DETECTION LIMITS


| | | |
|-----------------------------------|---------|---------|
| | LA/OA | C |
| Loading on primary filter (s/mm2) | 6.8E+01 | 6.8E+01 |
| Sensitivity (s/cc) | 2.2E-02 | 2.2E-02 |

****DO NOT USE AHERA UPLOAD SHEET****

Note: Currently, the spreadsheet will incorrectly exclude one particular class of matrix structure from the AHERA binned results summary -- a matrix for which a protruding bundle has an aspect ratio less than 5:1, but the underlying fibrils meet the aspect ratio requirement. This is a limitation of the TEM spreadsheet. If this type of structure is encountered, the analyst should note this issue in the analysis comments.

| Mineral Class | S<5u | S>5u | Q_ASBL | ASBL | Q_ASBC | ASBC |
|---------------|------|------|--------|------|--------|------|
| | | | | | | |

| Sensitivity |
|-------------|
| |



LIBBY**TEM Asbestos Structure Count -- ISO 10312****SAMPLE ID**

Status Analyzed
 EPA Sample Number FC-00007
 QA Type Not QA
 Lab Sample Number 143576-176058
 Sample Type Air
 Category Field
 Prep Indirect - Ashed
 Counting Rules 10312

PARAMETERS

Effective filter area 346.0 mm²
 Indirect factor 1.00E-01
 Number of Grid Openings (amphibole) 12
 Number of Grid Openings (chrysotile) 12
 Grid opening area 0.0110 mm²
 Volume (L) or Area (cm²) 1200 L
 Sensitivity (amphibole) 2.18E-02 s/cc
 Sensitivity (chrysotile) 2.18E-02 s/cc

Recording Rules:

| Min. AR | Min length (um) | Min width (um) |
|---------|-----------------|----------------|
| ≥ 3:1 | 0.5 | 0 |

Stopping Rules:

| Target S | Max GOs | Max N |
|----------|---------|-------|
| 0.001 | 100 | 50 |

COUNTS (based on countable structures only)

| Bin | LA | OA | C | PCME(all) | PCME(asb) |
|-------|----|----|---|-----------|-----------|
| a | 10 | 0 | 0 | | |
| b | 0 | 0 | 0 | | |
| c | 11 | 0 | 0 | | |
| d | 19 | 0 | 2 | | |
| e | 6 | 0 | 0 | | |
| f | 5 | 0 | 0 | | |
| Total | 51 | 0 | 2 | 17 | 17 |

Check OK OK OK

Grand total 53 OK

CONCENTRATION (s/cc)

| Bin | LA | OA | C | PCME(all) | PCME(asb) |
|-------|----------|-----|----------|-----------|-----------|
| a | 2.18E-01 | <DL | <DL | | |
| b | <DL | <DL | <DL | | |
| c | 2.40E-01 | <DL | <DL | | |
| d | 4.15E-01 | <DL | 4.37E-02 | | |
| e | 1.31E-01 | <DL | <DL | | |
| f | 1.09E-01 | <DL | <DL | | |
| Total | 1.11E+00 | <DL | 4.37E-02 | 3.71E-01 | 3.71E-01 |

| Type | Bin | Length | Width | Aspect ratio |
|---------------------------|-----|-------------|-------|--------------|
| LA = Libby-type amphibole | a | | | <5 |
| OA = Other amphibole | b | <.5 | | >= 5 |
| C = Chrysotile | c | | >.5 | >= 5 |
| | d | >=.5 to < 5 | <=.5 | >= 5 |
| | e | 5 to 10 | <=.5 | >= 5 |
| | f | >10 | <=.5 | >= 5 |

PCME: Length > 5 um, Width >= 0.25 um, Aspect Ratio >= 3:1

| Mineral Class | A_count | B_count | C_count | D_count | E_count | F_count | G_count |
|---------------|---------|---------|---------|---------|---------|---------|---------|
| LA | 10 | 0 | 11 | 19 | 6 | 5 | 51 |
| OA | 0 | 0 | 0 | 0 | 0 | 0 | 0 |
| C | 0 | 0 | 0 | 2 | 0 | 0 | 2 |

| Sensitivity | A_conc | B_conc | C_conc | D_conc | E_conc | F_conc |
|-------------|----------|--------|----------|----------|----------|----------|
| 2.18E-02 | 2.18E-01 | <DL | 2.40E-01 | 4.15E-01 | 1.31E-01 | 1.09E-01 |
| 2.18E-02 | <DL | <DL | <DL | <DL | <DL | <DL |
| 2.18E-02 | <DL | <DL | <DL | 4.37E-02 | <DL | <DL |

| G_conc |
|----------|
| 1.11E+00 |
| <DL |
| 4.37E-02 |