

REPORT

DEVELOPMENT AND VALIDATION OF AN ANALYTICAL METHOD FOR



**NOTOX Project 338805
NOTOX Substance 111834/B**

STATEMENT OF GLP COMPLIANCE

NOTOX B.V., 's-Hertogenbosch, The Netherlands

The study described in this report has been correctly reported and was conducted in compliance with the most recent edition of:

The OECD Principles of Good Laboratory Practice

which are essentially in conformity with:

The United States Food and Drug Administration. Title 21 Code of Federal Regulations Part 58.

The United States Environmental Protection Agency (FIFRA). Title 40 Code of Federal Regulations Part 160.


The United States Environmental Protection Agency (TSCA). Title 40 Code of Federal Regulations Part 792.

CONFIDENTIALITY STATEMENT

This report contains the unpublished results of research sponsored by [REDACTED] Chemicals B.V. Reproduction, issue or disclosure to third parties in any form is not permitted without prior written authorisation from the sponsor.

Study Director

[REDACTED]




Date: 18 July 2002

Management

[REDACTED]
Section Head

Analytical & Physical Chemistry



Date: July 18, 2002

QUALITY ASSURANCE STATEMENT

NOTOX B.V., 's-Hertogenbosch, The Netherlands

This report was audited by the NOTOX Quality Assurance Unit to ensure that the methods and results accurately reflect the raw data.

The dates of Quality Assurance inspections and audits are given below.
During the on-site inspections procedures applicable to this type of study were inspected.

DATES OF QAU INSPECTIONS/AUDITS

REPORTING DATES

on-site inspection (s)

18 February – 01 March 2002
(Process, physical chemistry)

04 March 2002

13 May – 31 May 2002
(Process, physical chemistry)

04 June 2002

protocol inspection (s)

06 November 2001 (Study)

06 November 2001

report audit (s)

13 June -18 June 2002 (Study)

18 June 2002

Head of Quality Assurance




Date: 19-7-02

SUMMARY

Analytical method

A high performance liquid chromatographic (HPLC) method for quantitative analysis of the test substance [REDACTED] was developed. A Zorbax RX-C18 column was used and a spectrophotometric detector set to read the absorbance at 260 nm for MIPK and 220 nm for all other components. During each run, the initial mobile phase of 46/54 (v/v) acetonitrile/Milli-Q water was changed gradually to 100 % acetonitrile by a gradient program.

Standard solutions of the test substance were prepared in acetonitrile at exactly known concentrations between 10950 mg/l and 23750 mg/l. In order to obtain suitable concentrations for the validation tests, these solutions were further diluted with 50/50 (v/v) acetonitrile/Milli-Q water.

Several small test substance peaks and one large test substance peak were observed in chromatograms of test substance solutions. The area of the small peaks was used as test substance response in calculations during the validation tests. The large peak corresponds to Dimethylphthalate which was used for preparation of this formulation. This peak was not used during method validation and not reported here since it cannot be analysed in the same run together with the smaller peaks due to the large difference in sensitivity. Peak identification was based on information supplied by the sponsor.

The HPLC method was validated for repeatability of injections, stability of the chromatographic system, stability of standard solutions, linearity and limit of detection.

Repeatability of injections based on [REDACTED]

The coefficient of variation of the response was determined to be 3.0 %, 1.6 % and 0.2 % for 10 replicates of a single 106 mg/l solution, a single 212 mg/l solution and a single 10575 mg/l solution, respectively. Therefore it was concluded that the chromatographic and integration conditions are acceptable for further analysis in the concentration range 212 - 10575 mg/l. Due to a relatively large coefficient of variation, it was necessary to pay extra attention to integration conditions at concentrations below 212 mg/l.

Repeatability of injections based on [REDACTED]

The coefficient of variation of the response was determined to be 0.2 %, 0.5 % and 0.2 % for 10 replicates of a single 106 mg/l solution, a single 212 mg/l solution and a single 10575 mg/l solution, respectively. Therefore it was concluded that the chromatographic and integration conditions are acceptable for further analysis in the concentration range 106 - 10575 mg/l.

Repeatability of injections based on [REDACTED]

The coefficient of variation of the response was determined to be 6.5 %, 3.9 % and 0.3 % for 10 replicates of a single 106 mg/l solution, a single 212 mg/l solution and a single 10575 mg/l solution, respectively. Due to a relatively large coefficient of variation, it was necessary to pay extra attention to integration conditions below 10575 mg/l.

Repeatability of injections based on [REDACTED] peak 1

The coefficient of variation of the response was determined to be 1.2 %, 0.7 % and 0.1 % for 10 replicates of a single 106 mg/l solution, a single 212 mg/l solution and a single 10575 mg/l solution, respectively. Therefore it was concluded that the chromatographic and integration conditions are acceptable for further analysis in the concentration range 106 - 10575 mg/l.

Repeatability of injections based on [REDACTED] peak 2

The coefficient of variation of the response was determined to be 0.8 %, 0.5 % and 0.1 % for 10 replicates of a single 106 mg/l solution, a single 212 mg/l solution and a single 10575 mg/l solution, respectively. Therefore it was concluded that the chromatographic and integration conditions are acceptable for further analysis in the concentration range 106 - 10575 mg/l.

Stability of the chromatographic system based on all components

The chromatographic system was stable over at least a 22.8-hour time interval at a test substance concentration of 10575 mg/l.

Stability of standard solutions based on all components

Standard solutions (11350 mg/l and 21460 mg/l) of the test substance in acetonitrile were stable for 24 days when stored at room temperature in the dark.

Linearity based on [REDACTED]

A linear relationship between response and test substance concentration was found over a concentration range of 21.2 – 10575 mg/l if a (1/concentration) weighting factor was applied ($n = 9$, $r = 0.9998$). Unweighted regression resulted in a deviation of more than 10% of the calibration points from the calculated line at concentrations below 212 mg/l.

Linearity based on [REDACTED] 4

A linear relationship between response and test substance concentration was found over a concentration range of 10.6 – 10575 mg/l if a (1/concentration) weighting factor was applied ($n = 10$, $r = 0.9998$). Unweighted regression resulted in a deviation of more than 10% of the calibration points from the calculated line at concentrations below 212 mg/l.

Linearity based on [REDACTED]

A linear relationship between response and test substance concentration was found over a concentration range of 52.9 – 10575 mg/l if a (1/concentration) weighting factor was applied ($n = 8$, $r = 0.9997$). Unweighted regression resulted in a deviation of more than 10% of the calibration points from the calculated line at concentrations below 212 mg/l.

Linearity based on [REDACTED] peak 1

A linear relationship between response and test substance concentration was found over a concentration range of 10.6 – 10575 mg/l if a (1/concentration) weighting factor was applied ($n = 10$, $r = 0.99994$). Unweighted regression resulted in a deviation of more than 10% of the calibration points from the calculated line at concentrations below 212 mg/l.

Linearity based on [REDACTED] peak 2

A linear relationship between response and test substance concentration was found over a concentration range of 10.6 – 10575 mg/l if a (1/concentration) weighting factor was applied ($n = 10$, $r = 0.99993$). Unweighted regression resulted in a deviation of more than 10% of the calibration points from the calculated line at concentrations below 212 mg/l.

Limit of detection based on [REDACTED] e

The limit of detection for the test substance was determined to be 0.8 mg/l at an injection volume of 100 μ l.

Limit of detection based on [REDACTED]

The limit of detection for the test substance was determined to be 2 mg/l at an injection volume of 100 μ l.

Limit of detection based on [REDACTED]

The limit of detection for the test substance was determined to be 46 mg/l at an injection volume of 100 μ l.

Limit of detection based on [REDACTED] peak 1

The limit of detection for the test substance was determined to be 11 mg/l at an injection volume of 100 μ l.

Limit of detection based on [REDACTED] peak 2

The limit of detection for the test substance was determined to be 10 mg/l at an injection volume of 100 μ l.

PREFACE

Sponsor

3

h

Study Monitor

D

SHERA, Regulatory Affairs

Testing Facility

NOTOX B.V.
Hambakenwetering 7
5231 DD 's-Hertogenbosch
The Netherlands

Study Director

Study plan

Start: 04 December 2001

Completed: 17 May 2002

TEST SUBSTANCE

Identification

Chemical name

CAS RN

Description

Clear colourless liquid

Batch

1510-14

Purity

See Certificate of Analysis

Test substance storage

In refrigerator in the dark

Stability under storage conditions

Stable

Expiry date

01 January 2003

The sponsor is responsible for all test substance data unless determined by NOTOX.

Note: **Don't heat up the test substance above 50°C**

PURPOSE AND PRINCIPLE

The purpose of the study was to implement and validate a method for the quantitative analysis of [REDACTED]. A high performance liquid chromatographic (HPLC) method was developed. The method was validated for repeatability of injections, stability of the chromatographic system, stability of standard solutions, linearity and limit of detection.

ARCHIVING

NOTOX B.V. will archive the following data for at least 10 years:
protocol, report, test substance reference sample and raw data.

REAGENTS

Milli-Q water	Tap water purified by reversed osmosis and subsequently passed over activated carbon and ion-exchange cartridges; Millipore, Bedford, MA, USA
Acetonitrile	HPLC-grade, Labscan, Dublin, Ireland

ANALYTICAL METHOD

A high performance liquid chromatographic (HPLC) method for quantitative analysis of the test substance was implemented by NOTOX B.V. The conditions are described below:

Column	Zorbax RX-C18, 250 x 4.6 mm; $d_p=5\ \mu\text{m}$ (Chrompack, Middelburg, the Netherlands)		
Mobile phase A	Acetonitrile		
Mobile phase B	Milli-Q water		
Gradient program	Time (min.)	%A	%B
	0	46	54
	5	46	54
	10	100	0
	13	100	0
	14	46	54
	19	46	54
Flow	2 ml/min		
Detection wavelength	220 nm and 260 nm		
Injection volume	100 μl		

PREPARATION OF SOLUTIONS

Standard solutions of the test substance were prepared in acetonitrile at exactly known concentrations between 10950 mg/l and 23750 mg/l. Suitable concentrations for the validation tests were obtained by dilution of these solutions with 50/50 (v/v) acetonitrile/Milli-Q water.

VALIDATION OF THE ANALYTICAL METHOD

The HPLC method was validated for repeatability of injections, stability of the chromatographic system, stability of standard solutions, linearity and limit of detection.

Repeatability of injections

Each of three solutions (106 mg/l, 212 mg/l and 10575 mg/l) was injected in tenfold into the HPLC system and the responses were recorded. The coefficient of variation of the response was calculated for each concentration.

Stability of the chromatographic system

A 10575 mg/l solution was injected (in duplicate) three times in a 22.8-hour time interval. The maximum deviation of the response was calculated.

Stability of standard solutions

Two standard solutions (11350 mg/l and 21460 mg/l) in acetonitrile were stored at room temperature in the dark and measured 24 days after preparation together with two freshly prepared standard solutions (10950 mg/l and 23750 mg/l). Prior to measurement, each of the standard solutions was diluted to an exactly known concentration of approximately 1100 mg/l. Each solution was injected in duplicate. The maximum deviation of the response factor was calculated.

Linearity

From one standard solution (21150 mg/l), ten dilutions were prepared. This resulted in a concentration range of 10.6-10575 mg/l. Each of the solutions was injected in duplicate. Responses were plotted against concentrations. A linear regression program was used to calculate the regression line from the responses and concentrations.

Limit of detection

A 10.6 mg/l and a 52.9 mg/l solution were injected in duplicate. In each chromatogram, the height (expressed in μ AU) of the test substance peak was measured as well as the noise level of the system (μ AU). The limit of detection was calculated from the mean peak height and the mean noise level.

DATA HANDLING

Response	Peak area of the test substance [units]
Response factor	Response/concentration [units* l * 10^3 /g]
Mean	$\bar{x} = \frac{1}{n} \sum_{i=1}^n x_i$ <p>where: x_i : measured value n : number of measurements</p>
Standard deviation	$s_{n-1} = \sqrt{\frac{\sum_{i=1}^n (x_i - \bar{x})^2}{(n-1)}}$
Coefficient of variation	(standard deviation/mean) * 100%
Maximum deviation	<p>[(highest - lowest)/mean] * 100%</p> <p>where 'mean' is the mean value of the highest and the lowest value.</p>
Linearity	<p>A linear regression program was used to calculate the regression line from the responses and concentrations. Linear regression analysis was performed using the least squares method. A (1/concentration) weighting factor was used.</p> <p>Regression line: $Y = aX + b$ where: Y : response X : concentration a : slope b : intercept</p>
Correlation coefficient	During regression analysis, the correlation coefficient (r) was calculated.
Limit of detection	<p>The limit of detection is defined as the lowest concentration of test substance that can be distinguished from instrumental noise using the analytical method described.</p> <p>Limit of detection = ((3 * noise level)/ signal) * conc.</p> <p>where: noise level (N) : height of the noise [μAU] signal (S) : height of the test substance peak [μAU] conc. : concentration of test substance [mg/l]</p>

RESULTS

A high performance liquid chromatographic (HPLC) method was developed for the quantitative analysis of [REDACTED]. The HPLC method was validated for repeatability of injections, stability of the chromatographic system, stability of standard solutions, linearity and limit of detection.

Several small test substance peaks and one large test substance peak were observed in chromatograms of test substance solutions. The area of all peaks was used as test substance response in calculations during the validation tests. Peak identification was based on information supplied by the sponsor. The large peak at $t = 3.11$ minutes corresponds to Dimethylphthalate which was used for preparation of this formulation. This peak was not used during method validation and not reported here since it cannot be analysed in the same run together with the smaller peaks due to the large difference in sensitivity.

The calculations for the validation tests were performed using not-rounded concentrations and responses. Therefore, some differences might be observed when calculating the statistical parameters using the values as mentioned in the tables.

Repeatability of injections

The results are summarised in Tables 1 - 5.

Table 1: Repeatability of injections based on Hydrogen peroxide.

Concentration [mg/l]	Response ¹ [units]	Standard deviation [units]	Coefficient of variation [%]
106	29780	893	3.0
212	45851	747	1.6
10575	1752466	3540	0.2

¹ Mean of 10 replicates of a single solution.

From these results, it was concluded that the chromatographic and integration conditions are acceptable for further analysis in the concentration range 106-10575 mg/l. Due to a relatively large coefficient of variation, it was necessary to pay extra attention to integration conditions at concentrations below 212 mg/l.

Table 2: Repeatability of injections based on [REDACTED]

Concentration [mg/l]	Response ¹ [units]	Standard deviation [units]	Coefficient of variation [%]
106	45640	111	0.2
212	89367	459	0.5
10575	4345925	7132	0.2

¹ Mean of 10 replicates of a single solution.

From these results, it was concluded that the chromatographic and integration conditions are acceptable for further analysis in the concentration range 106-10575 mg/l.

Table 3: Repeatability of injections based on [REDACTED].

Concentration [mg/l]	Response ¹ [units]	Standard deviation [units]	Coefficient of variation [%]
106	923	60	6.5
212	1940	77	3.9
10575	88680	309	0.3

¹ Mean of 10 replicates of a single solution.

Due to a relatively large coefficient of variation, it is necessary to pay extra attention to integration conditions below 10575 mg/l.

From these results, it was concluded that the chromatographic and integration conditions are acceptable for further analysis in the concentration range 106-10575 mg/l.

Table 4: Repeatability of injections based on [REDACTED] peak 1.

Concentration [mg/l]	Response ¹ [units]	Standard deviation [units]	Coefficient of variation [%]
106	10636	127	1.2
212	20600	153	0.7
10575	1015428	1446	0.1

¹ Mean of 10 replicates of a single solution.

From these results, it was concluded that the chromatographic and integration conditions are acceptable for further analysis in the concentration range 106-10575 mg/l.

Table 4: Repeatability of injections based on [REDACTED] peak 2.

Concentration [mg/l]	Response ¹ [units]	Standard deviation [units]	Coefficient of variation [%]
106	10077	79	0.8
212	19738	94	0.5
10575	984196	1355	0.1

¹ Mean of 10 replicates of a single solution.

From these results, it was concluded that the chromatographic and integration conditions are acceptable for further analysis in the concentration range 106-10575 mg/l.

Stability of the chromatographic system

The results are summarised in Tables 6 - 10.

Table 6: Stability of the chromatographic system based on [REDACTED]

Elapsed time [hours]	Response ¹ [units]	Maximum deviation [%]
0.0	1800553	3.0
0.4	1790596	
11.6	1746782	
11.9	1748878	
22.5	1753119	
22.8	1754237	

¹ For a 10575 mg/l solution.

Table 7: Stability of the chromatographic system based on [REDACTED]

Elapsed time [hours]	Response ¹ [units]	Maximum deviation [%]
0.0	4475698	3.3
0.4	4462955	
11.6	4331736	
11.9	4339414	
22.5	4363974	
22.8	4364589	

¹ For a 10575 mg/l solution.

Table 8: Stability of the chromatographic system based on [REDACTED]

Elapsed time [hours]	Response ¹ [units]	Maximum deviation [%]
0.0	90618	3.6
0.4	90126	
11.6	88974	
11.9	88360	
22.5	87419	
22.8	87422	

¹ For a 10575 mg/l solution.

Table 9: Stability of the chromatographic system based on [REDACTED] peak 1.

Elapsed time [hours]	Response ¹ [units]	Maximum deviation [%]
0.0	1048761	3.6
0.4	1050283	
11.6	1012958	
11.9	1015333	
22.5	1021622	
22.8	1022251	

¹ For a 10575 mg/l solution.

Table 10: Stability of the chromatographic system based on [REDACTED] peak 2.

Elapsed time [hours]	Response ¹ [units]	Maximum deviation [%]
0.0	1014356	3.4
0.4	1015497	
11.6	982026	
11.9	984518	
22.5	989923	
22.8	990744	

¹ For a 10575 mg/l solution.

From these results, it was concluded that the chromatographic system was stable over at least a 22.8-hour time interval at the concentration tested (10575 mg/l) for all measured components.

Stability of standard solutions

The results are summarised in Tables 11 - 15.

Table 11: Stability of standard solutions in acetonitrile based on [REDACTED]

Preparation date	Measurement date	Concentration ¹ [mg/l]	Response factor [units* 10^3 /g]	Maximum deviation [%]
23-4-02	17-5-02	11350	167.67/167.27	6.4
23-4-02	17-5-02	21460	178.07/177.80	
29-4-02	17-5-02	12500	178.48/176.52	
29-4-02	17-5-02	21150 ²	159.00/184.06	
16-5-02	17-5-02	10950	167.94/168.08	
16-5-02	17-5-02	23750	177.75/177.95	

¹ Concentration of the undiluted standard solution.

² Due to a large variation in the duplicate measurements, this solution was not used in the calculation.

From these results, it was concluded that standard solutions of the test substance (11350 mg/l and 21460 mg/l) in acetonitrile are stable for at least 24 days when stored at room temperature in the dark based on Hydrogen peroxide.

Table 12: Stability of standard solutions in acetonitrile based on [REDACTED]

Preparation date	Measurement date	Concentration ¹ [mg/l]	Response factor [units* 10^3 /g]	Maximum deviation [%]
23-4-02	17-5-02	11350	452.46/453.36	2.8
23-4-02	17-5-02	21460	442.22/442.77	
29-4-02	17-5-02	12500	441.27/441.96	
29-4-02	17-5-02	21150	452.51/453.30	
16-5-02	17-5-02	10950	440.73/440.88	
16-5-02	17-5-02	23750	446.86/447.16	

¹ Concentration of the undiluted standard solution.

From these results, it was concluded that standard solutions of the test substance (11350 mg/l and 21460 mg/l) in acetonitrile are stable for at least 24 days when stored at room temperature in the dark for MIPKP-T4.

Table 13: Stability of standard solutions in acetonitrile based on [REDACTED].

Preparation date	Measurement date	Concentration ¹ [mg/l]	Response factor [units*10 ³ /g]	Maximum deviation [%]
23-4-02	17-5-02	11350	9.8678/9.9057	14.9
23-4-02	17-5-02	21460	9.7661/9.8285	
29-4-02	17-5-02	12500	8.8296/8.7984	
29-4-02	17-5-02	21150	9.1584/9.3475	
16-5-02	17-5-02	10950 ²	7.2895/7.2639	
16-5-02	17-5-02	23750	8.6185/8.5293	

¹ Concentration of the undiluted standard solution.

² Due to a large variation in the solutions prepared on 16-5-02, this solution was not used in the calculation.

From these results, it was concluded that standard solutions of the test substance (11350 mg/l and 21460 mg/l) in acetonitrile are not stable for at least 24 days when stored at room temperature in the dark for [REDACTED]. Based on the results from other components it was assumed that test substance solutions are also stable for [REDACTED].

Table 14: Stability of standard solutions in acetonitrile based on [REDACTED] peak 1.

Preparation date	Measurement date	Concentration ¹ [mg/l]	Response factor [units*10 ³ /g]	Maximum deviation [%]
23-4-02	17-5-02	11350	105.26/105.56	2.5
23-4-02	17-5-02	21460	102.92/103.17	
29-4-02	17-5-02	12500	102.92/103.01	
29-4-02	17-5-02	21150	105.01/105.19	
16-5-02	17-5-02	10950	103.47/103.40	
16-5-02	17-5-02	23750	104.76/104.67	

¹ Concentration of the undiluted standard solution.

From these results, it was concluded that standard solutions of the test substance (11350 mg/l and 21460 mg/l) in acetonitrile are stable for at least 24 days when stored at room temperature in the dark for [REDACTED] peak 1.

Table 15: Stability of standard solutions in acetonitrile based on [REDACTED] peak 2.

Preparation date	Measurement date	Concentration ¹ [mg/l]	Response factor [units* 10^3 /g]	Maximum deviation [%]
23-4-02	17-5-02	11350	101.43/101.79	2.2
23-4-02	17-5-02	21460	99.612/99.980	
29-4-02	17-5-02	12500	99.585/99.809	
29-4-02	17-5-02	21150	101.84/101.84	
16-5-02	17-5-02	10950	99.737/99.819	
16-5-02	17-5-02	23750	101.19/101.16	

¹ Concentration of the undiluted standard solution.

From these results, it was concluded that standard solutions of the test substance (11350 mg/l and 21460 mg/l) in acetonitrile are stable for at least 24 days when stored at room temperature in the dark for [REDACTED] peak 2.

Linearity

The results are summarised in Tables 16 - 20, the corresponding regression line is shown in Figures 5 - 9.

Table 16: Linearity based on [REDACTED]

Concentration [mg/l]	Response ¹ [units]
10.6	17537 ² /17656 ²
21.2	18486/18552
52.9	22839 ² /24592
106	31231/30975
212	50094/49567
529	108058/106301
1058	190400/190454
2115	374666/375360
4230	732589/733806
10575	1746919/1749586
Slope	1.66×10^2
Intercept with Y-axis	1.51×10^4
Weighting factor	1/concentration
r	0.9998

¹ Duplicate measurements.

² Not used in calculations due to a large deviation from the calibration curve.

From these results, it was concluded that there is a linear relationship between response and concentration in the concentration range of 21 - 10575 mg/l if a (1/concentration) weighting factor is used for [REDACTED]. Unweighted regression resulted in a deviation of more than 10% of the calibration points from the calculated line at concentrations below 212 mg/l.

Table 17: Linearity based on [REDACTED]

Concentration [mg/l]	Response ¹ [units]
10.6	4884/4915
21.2	9358/9230
52.9	22804/22720
106	44655/44540
212	90295/90078
529	229192/228603
1058	447117/447566
2115	910092/910385
4230	1805880/1805044
10575	4343075/4350045
Slope Intercept with Y-axis Weighting factor r	4.18×10^2 7.04×10^2 1/concentration 0.9998

¹ Duplicate measurements.

From these results, it was concluded that there is a linear relationship between response and concentration in the concentration range of 11 - 10575 mg/l if a (1/concentration) weighting factor is used for [REDACTED]
[REDACTED] points from the calculated line at concentrations below 212 mg/l.

Table 18: Linearity based on [REDACTED]

Concentration [mg/l]	Response ¹ [units]
52.9	634/607
106	980/995 ²
212	1987/1996
529	4636/4708
1058	9276/9453
2115	18816/18757
4230	36792/36691
10575	88469/87803
Slope Intercept with Y-axis Weighting factor r	8.47×10^0 1.59×10^2 1/concentration 0.9997

¹ Duplicate measurements.

² Not used in calculations due to a large deviation from the calibration curve.

From these results, it was concluded that there is a linear relationship between response and concentration in the concentration range of 53 - 10575 mg/l if a (1/concentration) weighting factor is used for [REDACTED]. Unweighted regression resulted in a deviation of more than 10% of the calibration points from the calculated line at concentrations below 212 mg/l.

Table 19: Linearity based on [REDACTED] peak 1.

Concentration [mg/l]	Response ¹ [units]
10.6	1069/1072
21.2	1942/1818 ²
52.9	5160/5324
106	10424/10334
212	20522/20631
529	52143/52170
1058	103362/103427
2115	209574/209617
4230	414365/414508
10575	1015995/1019146
Slope	9.71×10^1
Intercept with Y-axis	5.82×10^1
Weighting factor	1/concentration
r	0.99994

¹ Duplicate measurements.

² Not used in calculations due to a large deviation from the calibration curve.

From these results, it was concluded that there is a linear relationship between response and concentration in the concentration range of 11 - 10575 mg/l if a (1/concentration) weighting factor is used for [REDACTED] peak 1. Unweighted regression resulted in a deviation of more than 10% of the calibration points from the calculated line at concentrations below 212 mg/l.

Table 20: Linearity based on [REDACTED] peak 2.

Concentration [mg/l]	Response ¹ [units]
10.6	959/960
21.2	1860/1768
52.9	4961/4880
106	10042/9844
212	19696/19768
529	50519/50530
1058	100478/100496
2115	203514/203672
4230	401266/401417
10575	985104/988086
Slope	9.42×10^1
Intercept with Y-axis	-5.46×10^1
Weighting factor	1/concentration
r	0.99993

¹ Duplicate measurements.

From these results, it was concluded that there is a linear relationship between response and concentration in the concentration range of 10.6 - 10575 mg/l if a (1/concentration) weighting factor is used for [REDACTED] peak 2. Unweighted regression resulted in a deviation of more than 10% of the calibration points from the calculated line at concentrations below 212 mg/l.

Limit of detection based on [REDACTED]

From two chromatograms of a 10.6 mg/l solution of the test substance, the mean noise level (N) was determined to be 0.06 mAU. From the same chromatograms, the test substance signal (S), i.e. the mean height of the test substance peak, was determined to be 2.39 mAU. Using these values, the limit of detection (S/N=3) was calculated to be 0.8 mg/l at an injection volume of 100 µl.

Limit of detection based on [REDACTED]

From two chromatograms of a 10.6 mg/l solution of the test substance, the mean noise level (N) was determined to be 0.06 mAU. From the same chromatograms, the test substance signal (S), i.e. the mean height of the test substance peak, was determined to be 0.93 mAU. Using these values, the limit of detection (S/N=3) was calculated to be 2 mg/l at an injection volume of 100 µl.

Limit of detection based on [REDACTED]

From two chromatograms of a 52.9 mg/l solution of the test substance, the mean noise level (N) was determined to be 0.03 mAU. From the same chromatograms, the test substance signal (S), i.e. the mean height of the test substance peak, was determined to be 0.103 mAU. Using these values, the limit of detection (S/N=3) was calculated to be 46 mg/l at an injection volume of 100 μ l.

Limit of detection based on [REDACTED] peak 1

From two chromatograms of a 10.6 mg/l solution of the test substance, the mean noise level (N) was determined to be 0.06 mAU. From the same chromatograms, the test substance signal (S), i.e. the mean height of the test substance peak, was determined to be 0.18 mAU. Using these values, the limit of detection (S/N=3) was calculated to be 11 mg/l at an injection volume of 100 μ l.

Limit of detection based on [REDACTED] peak 2

From two chromatograms of a 10.6 mg/l solution of the test substance, the mean noise level (N) was determined to be 0.06 mAU. From the same chromatograms, the test substance signal (S), i.e. the mean height of the test substance peak, was determined to be 0.185 mAU. Using these values, the limit of detection (S/N=3) was calculated to be 10 mg/l at an injection volume of 100 μ l.




Figure 1: HPLC chromatogram of a 10575 mg/l solution of XXXXXXXXXX 8 in 50/50 (v/v) acetonitrile/Milli-Q water at 220 nm [res.id. 1705].




Figure 2: HPLC chromatogram of a 10575 mg/l solution of XXXXXXXXXX in 50/50 (v/v) acetonitrile/Milli-Q water at 260 nm [res.id. 1704].




Figure 3: HPLC chromatogram of a blank solution (50/50 (v/v) acetonitrile/Milli-Q water) at 220 nm [res.id. 1703].




Figure 4: HPLC chromatogram of a blank solution (50/50 (v/v) acetonitrile/Milli-Q water) at 260 nm [res.id. 1702].




Figure 5: Regression line: response of [redacted] as a function of concentration [cal.curve id. 1847].




Figure 6: Regression line: response of [redacted] as a function of concentration [cal.curve id. 1688].




Figure 7: Regression line: response of [redacted] as a function of concentration [cal.curve id. 1687].




Figure 8: Regression line: response of [redacted] peak 1 as a function of concentration [cal.curve id. 1849].






Figure 9: Regression line: response of [redacted] peak 2 as a function of concentration [cal.curve id. 1690].




CERTIFICATE OF ANALYSIS

Certificate of AnalysisTNA-2001007
page 1 of 2

ICS-331

Product name : 
Chemical name : 
Batch number : 1510-14

Test results:


Method	Analysis of	Unit	Result *1
Jo/72.11,			
	Specification		
J20010792	D  e	% m/m	67.0 (± 1.0)
J20010792		% m/m	2.0 (± 0.3)
Amp/88.9	Water	% m/m	2.6 (± 0.3)
J20010792	Unidentified impurities	% m/m	0.5 (± 0.2)

*1 bracketed values are estimated 95% confidence intervals

File code : TNA-2001007

Analytical documentation : 20010792

Authorized by

Name : D 
Function : Section Head, Analytical Research Department
Date : October 25, 2001

Signature :



[REDACTED]

[REDACTED]

Certificate of Analysis

[REDACTED]

TNA-2001007
page 2 of 2

[REDACTED]

structure	% m/m
<div><div></div><div>(Type IV) IUPAC : [REDACTED]</div></div>	18.6
<div><div></div><div>(Type III) IUPAC : [REDACTED]</div></div>	7.9
<div><div></div><div>[REDACTED]e</div></div>	2.1

[REDACTED] 1.