

**U.S. Environmental Protection Agency  
Office of Research and Development**

**National Exposure Research Laboratory**

Research Triangle Park, North Carolina  
Athens, Georgia  
Cincinnati, Ohio  
Las Vegas, Nevada

**STANDARD OPERATING PROCEDURE (SOP)**

**Title: Analytical method for non-targeted and suspect screening in environmental and biological samples using Time of Flight Mass Spectrometry (TOFMS)**

**SOP Number: D-EMMD-PHCB-034-SOP-01**

**Effective Date: April 6, 2016**

**SOP was Developed:** ☒ **In-House** ☐ **Extramurally**

*Alternative Identification:*

**SOP Steward**

**Name: Mark Strynar**

**Signature / Date:**

**Approval**

**Name: Myriam Medina-Vera**

**Title: Branch Chief**

**Signature / Date:**

**Concurrence\***

**Name: Sania Tong-Argao**

**Title: Quality Assurance Manager (EMMD)**

**Signature / Date:**

\* To be used by Project Officers for review by the QAM or for Category 'A' projects.

**Disclaimer:** This SOP is an internal document intended for use in NERL laboratories for research purposes only. This is NOT an officially approved EPA Method. Any modification of or deviation from this SOP is documented according to Agency policy. Mention of brand names or vendors does not constitute an endorsement of products or services by the U.S. EPA.

## Analytical method for non-targeted and suspect screening in environmental and biological samples using Time of Flight Mass Spectrometry (TOFMS)

### Table of Contents

1.0	Scope and Application.....	2
2.0	Summary of Method.....	2
3.0	Definitions, Acronyms, and Abbreviations .....	3
4.0	Health and Safety .....	5
5.0	Cautions.....	5
5.1	Contamination of Samples.....	5
5.2	Sample Loss.....	5
6.0	Interferences .....	5
7.0	Personnel Qualifications/Responsibilities .....	6
8.0	Equipment and Supplies .....	6
9.0	Procedure.....	7
9.1	Cleaning of Labware and Surfaces .....	7
9.2	Overview of Sample Analysis .....	7
9.3	Placing Samples in Autosampler Vials .....	8
9.4	HPLC/TOFMS Analysis .....	8
9.4.1	HPLC/TOFMS Analysis .....	9
9.4.2	Tuning Solutions and Reference solutions for TOFMS Analysis.....	15
9.5	Batch order and equilibration .....	16
9.6	Standard Calibration Curve and Calculations.....	16
9.7	Nontargeted Screening and Suspect Screening.....	17
9.7.1	Molecular Feature Extraction .....	17
9.7.1.1	Method Explorer: Find Compounds-Find By Molecular Feature .....	18
9.7.1.2	Method Explorer: Find Compounds by Formula-Find By Formula Options.....	22
9.7.1.3	Method Explorer: Identify Compounds-Search Database .....	26
9.7.1.4	Method Explorer: Identify Compounds-Generate Formulas .....	30
9.8	Waste Management.....	32
9.8.1	Liquid Waste.....	32
9.8.2	Solid Waste.....	32
10.0	Data and Records Management .....	33
10.1	Laboratory Record Book .....	33
10.2	Location of LRB .....	33
11.0	Quality Control and Quality Assurance .....	33
12.0	References .....	33

## **Analytical method for non-targeted and suspect screening in environmental and biological samples using Time of Flight Mass Spectrometry (TOFMS)**

### **1.0 Scope and Application**

This standard operating procedure (SOP) describes a method for the non-targeted analysis (NTA) and suspect screening of chemicals in environmental and biological media using Time of Flight Mass Spectrometry (TOFMS). This method will be used to discover xenobiotic chemicals in various environmental and biological media screening against various in-house and OEM purchased data bases. There is much interest in the determination of the chemicals to get a more holistic idea of what we are exposed to in various media that are pertinent to human exposure studies (water, house-dust, soil, serum etc.). Classically this has been done through target analysis of a small number of chemicals of interest in select media. However, newer high-resolution mass spectrometry (HR-MS) techniques (such as TOFMS) allow for the rapid screening of analytes of interest based accurate mass measurements and spectral isotope matching. As the technique is kept consistent through the screening a various media, the result will be the ability to detect chemicals in different media, as well as co-occurrence of chemicals in specific media. This data will help to prioritize existing high-throughput toxicology assays (via NCCT) and modeling efforts (NERL) for chemicals of interest.

### **2.0 Summary of Method**

This method describes the LC/TOFMS analysis for the NTA and suspect screening of chemicals in various environmental and biological media. Refer to previously written SOPs for sample extraction and cleanup techniques. The following SOPs will be followed with minor alterations for NTA and suspect screening efforts. The SOPs include (<sup>1</sup>EMAB-112.0- Extraction and Analysis of Perfluorinated Compounds (PFCs) from Surface Waters by High Performance Liquid Chromatography (HPLC)-Tandem Mass Spectrometry (MS/MS); <sup>2</sup>MDAB-23.0- Standard Operating Procedure (SOP) For Extraction And High Performance Liquid Chromatography/Mass Spectrometry Analysis Of Perfluorinated Acids And Sulfonates From House Dust; and <sup>3</sup>MDAB 077.0- Extraction of Perfluorinated Compounds from Fish Liver Homogenate) for extraction and cleanup techniques. Additional information about the analytical equipment and software details can be found in the manuals (see References).

A portion of the sample extract is added to an auto-sampler vial along with 2 mM ammonium acetate buffer to match starting conditions (75% water: 25% methanol) of the HPLC gradient and is ready for LC/TOFMS analysis.

### 3.0 Definitions, Acronyms, and Abbreviations

- **Extract:** The sample extract that contains native target analytes and internal Standards.
- **Internal Standard (IS):** The fixed amount of compound that is added to each sample prior to its extraction. If quantitation is to be done the ratio of the detector signal of the native analyte to the detector signal of the IS is compared to ratios obtained from calibration curves where the IS level remains fixed and the native analyte levels vary. The IS is used to correct for minor sample-to-sample differences in extraction, matrix, purification, injection volume, chromatographic behavior, and MS ionization efficiency. All internal standards used in this procedure are mass-labeled, meaning that one or more atoms in the internal standard are replaced with a lesser abundant, higher molecular weight isotope, resulting in a mass difference between the analyte and mass labeled internal standard that allows the mass spectrometer to distinguish between the analyte and internal standard. In particular with NTA and suspect screening analysis, IS chemicals may be used to adjust for chromatographic drift from sample to sample, and to act as additional compounds to judge the mass accuracy of the instrumentations performance.
- **Double Blank:** A blank sample consisting of 100  $\mu$ L of MeOH and 300  $\mu$ L of 2 mM ammonium acetate buffer. Double blanks do not undergo complicated sample preparation. Adding the MeOH and buffer to an auto-sampler vial is the only sample preparation required. No internal standard or analytes are added to double blanks. This sample is used to assess solvent, and HPLC vial and cap contribution to analytical contamination.
- **Method Blank:** A blank sample that has undergone all sample preparation (i.e. addition of internal standard, alkaline digestion, and solid phase extraction (SPE)) as water, house-dust, soil or serum samples, but which began with DI water rather than an actual sample.
- **Matrix Blank:** A blank sample that has undergone all sample preparation (i.e. addition of internal standard, alkaline digestion, and SPE) as water, house-dust, soil or serum samples, if a blank matrix is available that is analyte free. In the event no blank matrix that is deemed acceptable, none will be used.
- **Quality Control (QC) Sample:** Depending on the media of choice this may include a standard reference material (SRM) from the National Institute of Standards and Technology (NIST) (house-dust, serum), or a prepared sample with a secondary standard source if available. QC samples undergo the same sample preparation (i.e. addition of internal standard, alkaline digestion, and SPE) as unknown samples. Measurements of quality control samples with each batch of samples gives the analyst a measure of whether or not sample preparation, extraction, and analysis has been performed properly.
- **Calibration standards:** In general calibration standards are not run in NTA and suspect screening efforts with respect to quantitation of target compounds in

environmental and biological media. However, with NTA and suspect screening efforts a series of calibrations compounds are used to mass calibrate the instrument daily before use and to auto-tune the TOFMS instrument. These consist of a series of known compounds in an OEM solution used to assure the mass accuracy of the instrument on a regular basis. In addition, a second set of known compounds (called reference compounds) are continually infused into the TOFMS for real-time mass correction. These are the compounds: purine [exact mass = 120.043596], HP0921 hexakis (1H,1H,3H-tetrafluoropropoxy) phosphazene [exact mass = 921.002522] and tetrahydroperfluorononanoic acid (THPNA) [exact mass = 391.0009]. Depending on the polarity of the instrument and the mobile phase modifiers used, different reference masses are seen. In addition, any known compound with an exact mass may be added that is not expected to be present in the samples being analyzed. Refer to Table 1 for additional references masses and forms used in this analysis.

Table 1: Reference masses for real-time mass correction in TOFMS analysis

Species	Positive Ion m/z	Negative Ion m/z
CF <sub>3</sub> (trifluoro acetic acid TFA) fragment)		68.995758
TFA anion		112.985587
purine	121.050873	119.036320
HP0921	922.009798	
HP0921 (formate adduct)		966.000725
HP0921 (acetate adduct)		980.016375
HP0921 (TFA adduct)		1033.988109
THPFNA		391.0009

- **Unknown Sample:** Each unknown sample is extracted according to the aforementioned SOPs and analyzed along with QA/QC, blank and replicate sample. This may consist of any environmental or biological media (i.e. house-dust, soil, water, serum...) and is compared to other samples with unknown chemical compositions. When available metadata (age, gender, location, upstream/downstream...) will be used to draw conclusions from comparisons to other unknown samples.
- **DI:** Deionized
- **SPE:** Solid phase extraction
- **LTL:** Laboratory team leader

- **LRB:** Laboratory record book
- **PPE:** Personal protective equipment

#### **4.0 Health and Safety**

Proper personal protective equipment (PPE) should be used when performing this procedure. PPE that should be worn throughout this procedure include nitrile gloves, safety glasses, and a lab coat. Solvents should be opened and handled in a properly functioning fume hood.

#### **5.0 Cautions**

##### **5.1 Contamination of Samples**

Care must be taken to avoid sample contamination. Throughout sample preparation, laboratory personnel must remove gloves that have been or are suspected to be contaminated and put on a clean pair of gloves. A clean pipette tip must be used for each sample. Solvents and DI water should be tested to regularly to ensure they are not significant sources of contamination. All glassware and surfaces used in sample preparation must be as clean as possible before usage. Blanks (solvent and method) must be run with each batch of samples to determine what contamination has occurred during sample preparation and analysis. Standards and samples must be capped when not in use to avoid evaporative loss of solvent (MeOH), which could lead to increased concentration of analytes. Standards and samples should not be stored in the same refrigerator or freezer.

##### **5.2 Sample Loss**

Loss of some analytes, especially lower carbon chain compounds may occur at high temperatures (greater than 60°C) and when samples are evaporated to dryness. Therefore, all sample preparation and analysis should occur at or below 60°C and samples should not be evaporated to dryness (Section 9.6 below). Another potential loss of chemicals is due to sorption to labware and precipitation due to cooling of samples in a refrigerator (4°C). As this in NTA and suspect screening, all unknown compounds cannot be assessed for losses. However, general laboratory practice based on past experience should be taken to help minimize losses. Assessment of losses of already known compounds in samples may be useful to help judge method performance. In addition the use of ISs (see section 3.0) may be beneficial in assessing losses as well.

#### **6.0 Interferences**

As HRMS techniques are monitoring for all m/z at the entire time of an analytical

run interferences are expected. This is why it is critical to run solvent and method blanks along with unknown and QA/QC samples. The presence of interferences in particular samples should be assessed by comparisons to blanks. This may be done in several ways. In particular the chemicals found in blank samples may be placed in an “exclusion list” to negate being found in unknown samples and occupy valuable analytical time. In addition, chemicals found in composite of solvent and process blanks can be combined into a personal compound database list (PCDL) with m/z to be sure to eliminate these chemicals based on both m/z and RT rather than simply m/z. This is important as a simple m/z exclusion list may exclude an important mass if it is similar to a chemical found in the blanks. An additional approach is to not exclude any chemicals and when done compare blank samples to unknown samples. Any peak in a sample found to have significantly more signal than in the blank (i.e. S/N > 3:1) could be deemed a real chemical.

## **7.0 Personnel Qualifications/Responsibilities**

Note that this SOP assumes a thorough working knowledge of basic laboratory skills, reagents, and instrumentation. This document is designed to guide a trained laboratory worker in the methodology discussed below and it is not intended to instruct individuals on the basic aspects of analytical chemistry. Performance of procedures described in this SOP requires previous laboratory experience and training. At least 2 years of previous laboratory experience is necessary to perform the procedures described in this SOP and a technical or 4-year degree in chemistry would be desirable. The laboratory worker who performs the sample extractions will be responsible for entering experimental information in the laboratory record books. The Laboratory Team Leader (LTL) will also oversee the sample preparation, extraction, and cleanup, to ensure all SOPs are followed by all project staff. However, this SOP only covers NTA and suspect screening analysis of environmental and biological samples.

## **8.0 Equipment and Supplies**

### **8.1 Laboratory Equipment**

- Fume hood
- Solvent cabinet
- Vortexer
- Nitrogen Generator or clean compressed Nitrogen gas supply (>99% purity)
- Agilent 1100 high pressure liquid chromatography system with autosampler (Agilent Technologies, Wilmington DE)

- HPLC column Agilent Eclipse Plus C18 (2.1 × 50 mm, 3.5 µm) (Agilent Technologies, Wilmington DE) or comparable
- Agilent TOF MSD model 6200 (Agilent Technologies, Wilmington DE)

## **8.2 Laboratory Supplies**

- Variable volume standard pipettors (10 – 100 µL, 100 – 1000 µL, and 500 – 5000 µL)
- Variable volume repeat pipettors (Eppendorf Repeater Plus)
- RNAase and DNAase free Pipette tips (Genessee Scientific, San Diego, CA)
- Sterile Eppendorf repeat pipettor tips (50 mL)
- Polypropylene autosampler vials and caps (Laboratory Supply Distributors Corp., Mt. Laurel, NJ)
- Ultra-high purity grade compressed nitrogen, (HoloX Co., Norcross, GA)
- Solvent bottles
- Permanent markers

## **8.3 Chemicals and Reagents**

- HPLC-grade methanol (MeOH), acetonitrile (ACN) and isopropyl alcohol from Burdick-Jackson (Muskegon, MI) or comparable
- Deionized water from a Barnstead EASYpure UV/UF compact reagent grade water system (Dubuque, IA)
- Ammonium acetate from Sigma-Aldrich Chemical (St Louis, MO)
- Ammonium formate from Sigma-Aldrich Chemical (St Louis, MO)
- Agilent 5 mM purine in acetonitrile:water Agilent part # 18720242
- Agilent 2.5 mM HP0921 in acetonitrile:water Agilent part # 18720241
- Agilent ESI-L Low concentration Tuning mix Agilent part # G1969-85000 (100 mL)

## **9.0 Procedure**

### **9.1 Cleaning of Labware and Surfaces**

- 9.1.1 Wash all glass equipment used in the procedure using a laboratory-grade dishwasher.
- 9.1.2 Rinse all glass equipment with methanol immediately before use.
- 9.1.3 Rinse beakers that will contain aqueous reagents once with deionized water.
- 9.1.4 Cover all laboratory bench surfaces used in the procedure with clean bench-top laboratory paper.

### **9.2 Overview of Sample Analysis**

Each sample analysis set consists of the following:



- Double blanks ( $n$  = dependent on number of samples to be run; at a minimum run before and after each set. In addition interspersed in between unknown samples not to exceed once every 10 unknown samples) mobile phase solvents at ratio equal to starting conditions of gradient.
- Method blank ( $n \geq 1$ ) (procedural blank)
- Matrix blanks ( $n \geq 1$ ) (if available; matrix blank should match unknown samples being analyzed)
- Calibration curve standards if applicable ( $n \geq 6$ ) (analytes of interest spiked into blank matrix if available; otherwise multiple extracted procedural blanks)
- QC samples ( $n$  = generally 10% of the number of unknown samples) (either QCs made from of a purchased SRM or spiked in to a blank matrix and run prior to assay to generate target statistics then with sequential analytical batches to assure quality control)
- Unknown samples (typically 1 to 50 samples, each consisting of extracts from environmental or biological samples)
- Replicates ( $n$  = generally 10% of the number of unknown samples)

All method and matrix blanks, QCs, calibration standards, replicate and unknown samples are subjected to the same sample preparation and analysis. The following procedure assumes that all method and matrix blanks, QCs, calibration standards, and unknown fish samples have already extracted and cleaned up via the SOPs mentioned in section 2.0.

### 9.3 Placing Samples in Autosampler Vials

The condensed organic solvent (generally methanol or acetonitrile) eluate is vortexed and placed in a polypropylene autosampler vial containing 2 mM ammonium acetate buffer. The ratio of solvent:buffer should approximate the starting conditions of the mobile phase of the gradient used. Each vial is then capped and vortexed to ensure mixing of the organic sample with the aqueous buffer. The bottom of each vial is checked for air bubbles, which can be removed by tapping on the vial. Samples are then placed in the HPLC autosampler and analyzed by LC/TOFMS.

### 9.4 HPLC/TOFMS Analysis

Samples were analyzed using an Agilent 1100 HPLC (Agilent Technologies, Palo Alto, CA) interfaced with an Agilent LC/MSD TOF model # G1969A (Agilent Technologies, Palo Alto, CA)). A 45-minute gradient HPLC run (Table 2) was used with a mobile phase components of MeOH and 2 mM formate or acetate buffer at a flow rate of 300  $\mu\text{L}/\text{min}$ . The HPLC column was an Agilent Eclipse Plus C18 (2.1  $\times$  50 mm, 3.5  $\mu\text{m}$ ) with an injection volume of 40  $\mu\text{L}$ . A comparable HPLC column may be substituted. Electrospray ionization was used in the mass spectrometer source, which was maintained at 325°C. Depending on

the polarity of the instrument and the mobile phase used, specific reference mass ions are monitored for during analytical runs due to continual infusion for mass correction. See Table 1 for reference ions  $m/z$ .

**Table 2.** The gradient program used on Agilent 1100 HPLC system for the NTA and suspect screening analysis of environmental samples should be kept consistent to allow for chromatographic separation of matrix interferences from analytes of interest and to allow for comparison between sample extracts from different media. However, this gradient may be altered to suit individual samples sets if needed.

Time (min)	Flow Rate (mL/min)	%A <sup>a</sup>	%B <sup>b</sup>
0	0.2	75	25
25	0.2	20	80
40	0.2	0	100
45	0.2	0	100
Post time (4 minutes)	0.2	75	25

<sup>a</sup>Mobile phase component A consisted of 2 mM ammonium format or acetate in deionized water

<sup>b</sup>Mobile phase component B consisted of methanol; acetonitrile may be used if desired for additional assay if needed.

#### 9.4.1 HPLC/TOFMS Analysis

Specific components of the TOFMS Methods are broken up into 4 tabs on the Mass Hunter software. The simplest way to replicate the methods employed are to replicate the setting in the software. The component tabs are:

- 1) uHiP ALS (the autosampler settings)
- 2) CapPump ( the HPLC pump settings)
- 3) TCC (the column compartment settings)
- 4) MS TOF ( the TOF mass spectrometer settings)

Within each method tab are several tabs for additional settings. This section is organized by component section, then by additional tabs as screen shots.

## uHiP ALS Tab

### 1. Setup Tab

Properties   <b>uHiP ALS</b>   CapPump   TCC   MS TOF	
Setup   Options   Carryover Reduction	
<p>Injection</p> <p><input type="radio"/> Standard Injection</p> <p><input checked="" type="radio"/> Injection with Needle Wash</p> <p><input type="radio"/> Use Injector Program      Total lines: -</p> <p>Injection Volume: <input type="text" value="40"/> <math>\mu</math>l</p>	<p>High Throughput</p> <p><input checked="" type="radio"/> Disable Overlapped Injection</p> <p><input type="radio"/> Start overlap when sample is flushed out</p> <p><input type="radio"/> Start overlap after <input type="text" value="0"/> minutes</p> <p><input type="checkbox"/> Automatic delay volume reduction</p>
<p>Needle Wash</p> <p><input type="radio"/> Wash the needle in vessel <input type="text" value="Vial 1"/> , and repeat <input type="text" value="1"/> times</p> <p><input checked="" type="radio"/> Wash needle in Flush Port for <input type="text" value="4"/> seconds</p>	

### 2. Options Tab

Properties   <b>uHiP ALS</b>   CapPump   TCC   MS TOF	
Setup   Options   Carryover Reduction	
<p>Stop Time</p> <p><input checked="" type="radio"/> As pump</p> <p><input type="radio"/> No limit</p> <p><input type="radio"/> <input type="text" value="0"/> min</p>	<p>Post Time</p> <p><input checked="" type="radio"/> Off</p> <p><input type="radio"/> <input type="text" value="0"/> min</p>
<p>Auxiliary</p> <p>Draw Speed: <input type="text" value="50"/> <math>\mu</math>l/min      Equilibration Time: <input type="text" value="0"/> sec</p> <p>Eject Speed: <input type="text" value="100"/> <math>\mu</math>l/min      Sample Flush-out factor: <input type="text" value="5"/></p> <p>Draw Position: <input type="text" value="0"/> mm      <input type="checkbox"/> Store temperature</p> <p><input type="checkbox"/> Vial/Well bottom sensing</p>	

### 3. Carryover Reduction Tab

Properties	<b>μHiP ALS</b>	CapPump	TCC	MS TOF
Setup	Options	<b>Carryover Reduction</b>		
Cleaning Wellplate Sampler				
Cleaning Injection Valve		<a href="#">Settings...</a>		
<input type="checkbox"/> Enable Rinse				
Rinse Draw Speed	100	μl/min	Rinse Volume (organic)	1.3 * volume
Rinse Eject Speed	100	μl/min	Rinse Volume (water)	2 * volume
(volume = syringe + loop cap.+ seat cap.)				
Please refer to online help when using HPLC Chip Carryover Reduction.				

## CapPump Tab

### 1. Setup Tab

Properties	μHiP ALS	<b>CapPump</b>	TCC	MS TOF
Setup	Timetable	Options		
<b>Mode</b> <input type="radio"/> Micro Flow <input checked="" type="radio"/> Normal Flow		<b>Stop Time</b> <input type="radio"/> No Limit <input checked="" type="radio"/> 45 min		
<b>Flow</b> Column Flow: 200 μl/min		<b>Post Time</b> <input type="radio"/> Off <input checked="" type="radio"/> 4 min		
<b>Calibrated as</b> <input checked="" type="radio"/> Uncalibrated (H2O-H2O) <input type="radio"/> H2O-DMF - A:Aqueous B:Dimethyl <input type="radio"/> Custom Calibration <a href="#">Setup</a>		<b>Pressure Limits</b> Min: 0 bar   Max: 400 bar		
		<b>Solvent A</b> 75.00 % 1: <input type="text"/> 2: 95:5 DI:MeOH 0.4 mL		
		<b>Solvent B</b> <input checked="" type="checkbox"/> 25 % 1: <input type="text"/> 2: 95:5 MeOH:DI 0.4 mL		
		<b>Fast Reconditioning</b> <input type="radio"/> On <input checked="" type="radio"/> Off		

## 2. Timetable Tab (same as Table 2 information above)

Properties   <u>μHiP ALS</u>   CapPump   TCC   MS TOF				
Setup   Timetable   Options				
	Time	B%	Flow	Max. Press.
1	0.00	25.0	200.000	
2	25.00	80.0	200.000	
3	40.00	100.0	200.000	
4	45.00	100.0	200.000	

## 3. Options Tab

Properties | μHiP ALS | CapPump | TCC | MS TOF

Setup | Timetable | Options

Maximum Flow Gradient

100 ml/min per minute

Store

☐ % A ☐ Flow

☐ % B ☐ Pressure

☐ Piston movements A

☐ Piston movements B

Primary Flow

Low solvent consumption Fast gradients

500 .. 800 μl/min

Minimum Stroke

Channel A: 20 μl

Channel B: 20 μl

Auto

Compressibility

Channel A: 50 \*10<sup>-6</sup>/bar

Channel B: 115 \*10<sup>-6</sup>/bar

Column Flow Ready Condition

☐ Always ready

☒ Within +/- 3 %

## TCC Tab

### 1.Setup Tab

Properties | μHiP ALS | CapPump | TCC | MS TOF

Setup | Timetable | Options

Temperature (Left)

☒ 35 °C

☐ Not controlled

Temperature (Right)

☐ 20 °C

☐ Not controlled

☒ Same as left

Stop Time

☒ As pump / injector

☐ No limit

☐ 0.1 min

Post Time

☒ Off

☐ 0 min

Column Switching Valve

Column 1

## 2.Options Tab

The screenshot shows the 'Options' tab within a software window. The window has tabs for 'Properties', 'μHiPALS', 'CapPump', 'TCC', and 'MS TOF'. The 'Options' tab is active, showing two sub-tabs: 'Setup' and 'Options'. The 'Options' sub-tab is selected, displaying two main sections: 'Store' and 'Enable Analysis'.

**Store Section:**

- ☒ Temperature (left)
- ☐ Temperature (right)

**Enable Analysis Section:**

- ☐ With any temperature
- ☒ When temperature is within setpoint

Below the 'When temperature is within setpoint' option, there is a field for the setpoint range:  $\pm 0.8$  °C.

## MS TOF Tab

### 1. General Tab

The screenshot shows the 'MS TOF' tab within a software window. The window has tabs for 'Properties', 'μHiPALS', 'CapPump', 'TCC', and 'MS TOF'. The 'MS TOF' tab is active, showing a 'General' sub-tab. The 'General' sub-tab displays various configuration options for the MS TOF system.

**General Section:**

- Ion Source:** Dual ESI
- Ion Polarity:** Positive
- Data Storage:** Both
- LC Stream:** MS

**Stop Time:**

- ☒ No Limit/As Pump
- ☐ Stop Time: 30 min

**Time Segment and Experiment #:**

- Time (min):** 0
- Expt #:** 1

**Ion Polarity (Seg):**

- ☒ Positive
- ☐ Negative
- ☐ Fast Polarity Switching

**Data Storage (Seg):**

- ☐ None
- ☒ Both
- ☐ Centroid
- ☐ Profile

**LC Stream (Seg):**

- ☒ MS
- ☐ Waste

**Plot and Centroid Data Storage Threshold:**

- Abs. threshold:** 200
- Rel. threshold (%):** 0.01

**Apply Now** button.

☐ Do not wait for setpoints (e.g. temperature) to equilibrate

## 2. Source Tab

Properties | **μHPLS** | CapPump | TCC | **MS TOF**

Ion Source: Dual ESI Ion Polarity: Positive Data Storage: Both LC Stream: MS

Stop Time: ☐ No Limit/As Pump ☐ Stop Time 30 min

Time Segment and Experiment #: Time (min) 0 Expt 1

Cycle Time 1 s

General | **Source** | Acquisition | Ref Mass | Chromatogram

Dual ESI (Seg): Gas Temp 350 °C 0 °C Drying Gas 10 l/min 0.0 l/min Nebulizer 30 psig 0 psig

MS TOF (Expt): Fragmentor 80 V Skimmer 65 V OCT 1 RF Vpp 250 V

Dual ESI (Expt): VCap 3500 V Capillary 0.000 uA Chamber 0.00 uA

## 3. Acquisition Tab

Properties | **μHPLS** | CapPump | TCC | **MS TOF**

Ion Source: Dual ESI Ion Polarity: Positive Data Storage: Both LC Stream: MS

Stop Time: ☒ No Limit/As Pump ☐ Stop Time 30 min

Time Segment and Experiment #: Time (min) 0 Expt 1

Cycle Time 1 s

General | Source | **Acquisition** | Ref Mass | Chromatogram

TOF Spectra: Mass Range Min Range 100 m/z Max Range 1700 m/z

Acquisition Rate/Time: Rate 1 spectra/s Time 1000 ms/spectrum Transients/spectrum 13682

## 4. Ref Mass Tab

Properties | **μHPLS** | CapPump | TCC | **MS TOF**

Ion Source: Dual ESI Ion Polarity: Positive Data Storage: Both LC Stream: MS

Stop Time: ☒ No Limit/As Pump ☐ Stop Time 30 min

Time Segment and Experiment #: Time (min) 0 Expt 1

Cycle Time 1 s

General | Source | Acquisition | **Ref Mass** | Chromatogram

Reference Mass Correction: ☒ Enable ☒ Use bottle A Apply Now Ref Nebulizer 15 psig

Auto Recalibration Reference Mass Parameters: Detection Window 100 ppm Minimum Height 1000 counts

Reference Masses: Reference Masses Table

On	M/Z
<input checked="" type="checkbox"/>	121.050873
<input type="checkbox"/>	149.02332
<input type="checkbox"/>	322.048121
<input checked="" type="checkbox"/>	922.009798
<input type="checkbox"/>	1221.990637
<input type="checkbox"/>	1521.971475
<input type="checkbox"/>	2421.91399

## 5. Chromatogram Tab

The screenshot shows the TOFMS software interface with the 'Chromatogram' tab selected. The interface includes several tabs at the top: Properties, uHPLC ALS, CapPump, TCC, MS TOF, General, Source, Acquisition, Ref Mass, and Chromatogram. The 'Chromatogram' tab is active, displaying a table with the following columns: Chromatogram, Label, Expt Type, Polarity Type, Offset, and Y-Range. The table contains one entry: TIC, TIC, MS, Both, 15, 10000000. Other settings visible include Ion Source: Dual ESI, Ion Polarity: Positive, Data Storage: Both, LC Stream: MS, Stop Time: No Limit/As Pump, and Cycle Time: 1 s.

### 9.4.2 Tuning Solutions and Reference solutions for TOFMS Analysis

Critical for any TOFMS or other HRMS analysis is the preparation of proper solutions for calibration of the instrument (tuning solution) and real-time mass correction during any analytical run (reference solutions). The tuning solutions and reference solutions are made up using Agilent stock solutions that can be purchased for an authorized vendor. See section 8.3 for further details.

#### Tuning Solution

- Agilent ESI-L Low concentration Tuning mix Agilent part # G1969-85000 (100 mL)

Positive Mode Dual ESI (use Undiluted)

Negative Mode Dual ESI (25 mL undiluted + 75 mL of 95:5 Acetonitrile:DI water)

#### Reference Solution

These solutions are required for real-time mass correction of all generated spectra. These are starting conditions for reference solutions. However, additional solution alteration may be required if the signal is too high or too low depending on the instrument's response and mode being used. It is critical that the reference solution  $m/z$  be able to be seen by the TOFMS without saturation of the detector. Generally, these values should be between 50,000 and 500,000 counts to be in this range.

- Agilent 5 mM purine in acetonitrile:water Agilent part # 18720242
- Agilent 2.5 mM HP0921 in acetonitrile:water Agilent part # 18720241
- \*\*\*additional reference compound added by RTP lab; Not Agilent reference



solution; 1000 ng/uL tetrahydro perfluorononanoic acid (THPNA)\*\*\*\*

**Positive Mode Dual ESI**

500 mL of Acetonitrile:DI water (90:10)  
1.5 mL of 5mM purine solution (see above)  
750 uL 2.5 mM HP0921solution (see above)

**Negative Mode Dual ESI**

1000 mL of Acetonitrile:DI water (90:10)  
300 uL of 5mM purine solution (see above)  
150 uL 2.5 mM HP0921solution (see above)  
100 uL of 1000 ng/uL solution of THPFNA

**9.5 Batch order and equilibration**

Prior to any analytical runs the HPLC/TOF MS is equilibrated for at least 15 minutes with the HPLC system running to waste and the TOF MS system in “On” mode. After equilibration a minimum of 2 double blanks are run prior to any analytical batch to ensure that mobile phase components and the instrument is sufficiently equilibrated. If contamination is determined, appropriate actions are taken to rectify the situation before analysis continues. Once the instrument is determined to be ready for analysis, samples are run in the order of double blanks, method blank, matrix blank, and calibration curve (if needed). After the highest calibration curve point a double blank is run prior to the next sample to ensure carry-over is not an issue. From here on out, QC samples, unknowns and duplicate samples are run in sets of 10 separated by a double blank until all samples have been run. After a double blank the entire calibration curve is then run again at the end of the analytical run (if needed), followed by multiple double blanks to prepare the system for the next analyst or assay.

**9.6 Standard Calibration Curve and Calculations**

In general calibration curves are not generally run with NTA and suspect screening analysis. However, at times it is necessary to calibrate for a particular compound of interest, or a series of compounds. This may be used as a semi-quantitative analysis as the analytes of interest generally far exceed in numbers the authentic standards available. In this instance it may be deemed necessary for a general calibration curve of multiple analytes be used for the calibration of analytes found via NTA or suspect screening. A calibration curve of 6-points should be sufficient to get adequate response over the linear response range of the analytes. However, since most analytes do not have authentic standards, conclusions drawn concerning quantitation should be used with caution.

In another application it is possible to use a one point calibration as a point estimate of compounds found via NTA or suspect screening. In this instance an

analyst could use a SRM or a mixed calibration solution with many compounds (100 or greater) to get a point estimate. In the mixed standard approach it is advisable for the standard to be run at two different concentrations (preferably serial dilution: 10:1) and not the response of the analyte being 10x between dilution. Again as before caution should be used in making conclusion about quantitative amounts, and should only be used when estimation of concentrations is needed.

## **9.7 Nontargeted Screening and Suspect Screening**

Generally after a set of samples is run on the HPLC TOFMS system the data processing needs to be accomplished offline through a series of procedures. These procedures are designed to look for features in a sample that match a chemical in a database (suspect screening), or to find chemicals in a sample for further structure elucidation (NTA). A feature is defined as (an accurate mass, a RT and an integrated peak area). The procedures that follow are the best practices to date that have been employed for each approach. However, each step in this procedure is open to further refinement and should be seen as a starting off point for post run analysis.

For NTA a molecular feature extraction must be run. For suspect screening one can skip the MFE step and use the Find by Formula algorithm, or run the suspect screening step (via identify compounds) against the chemicals found via MFE. If the MFE step is not used prior to performing suspect screening, the analysts will only ever find chemicals listed in the databases used. In this application, if that is acceptable suspect screening depends upon the database chosen. A culmination of the chemicals found from the various databases used (DSS-TOX, PCDLs and .CSV files) can be prepared in bulk.

The Generate Formula algorithm can be used to produce candidate formulas for a spectral peak or for compounds found via the MFE. The user may input elements of choice via the “allowed species” tab. Refer to section 9.7.1.4 for screenshots of this procedure. Once all of the parameters for formula generation are selected, a target spectra selecting the first spectral peak of the cluster via left mouse click and drag is chosen and the “play button” is selected to generate a formula for the spectral cluster. The first ion in the cluster must be the monoisotopic mass peak. Generally this is the largest peak. However, in some instances (high chlorinated or brominated compounds) third peak in the cluster is the largest.

### **9.7.1 Molecular Feature Extraction**

The intention of the MFE algorithm is to identify “features” in the total ion chromatogram (TIC) that require further attention. This MFE process has many variables that may be changed. This procedure is the currently employed set of variables that have been used. However, many of these variables may be changed

to user specific settings for a given project's application. For instance, one could tradeoff speed and detail by altering parameters such as "minimum signal intensity threshold" or "only viewing the top 100 features". It is dependent upon the analyst to determine if the parameters should be altered.

It is highly advisable that the data processing method be saved and dated in a similar fashion as the data files being processed. In this instance, the analyst can go back in the future and observe the parameters that were chosen for the post-hoc analysis.


Example:

Worklist: 4-18-16 Cape Fear water sample

Data Files: 4-18-16 Cape Fear water samples D.01 – D0x.

Data Analysis Method: 4-18-16 Cape Fear water data analysis

#### 9.7.1.1 Method Explorer: Find Compounds-Find By Molecular Feature



## Extraction Tab

Extraction	Ion Species	Charge State	Compound Filters	Mass Filters	Mass Defect	Results	Advanced
------------	-------------	--------------	------------------	--------------	-------------	---------	----------

Extraction algorithm

Target data type: Small molecules (chromatographic)

Input data range

☐ Restrict retention time to  minutes

☐ Restrict m/z to  m/z

Peak filters

☐ Use peaks with signal-to-noise  $\geq$   (Profile spectra only)

☒ Use peaks with height  $\geq$   counts (Profile and centroid spectra)

## Ion Species Tab

Extraction	Ion Species	Charge State	Compound Filters	Mass Filters	Mass Defect	Results	Advanced
------------	-------------	--------------	------------------	--------------	-------------	---------	----------

Allowed ion species

Positive ions	Negative ions	Neutral losses
<input checked="" type="checkbox"/> +H <input checked="" type="checkbox"/> +Na <input checked="" type="checkbox"/> +K <input checked="" type="checkbox"/> +NH <sub>4</sub>	<input checked="" type="checkbox"/> -H <input type="checkbox"/> +Cl <input type="checkbox"/> +Br <input checked="" type="checkbox"/> +HCOO <input type="checkbox"/> +CH <sub>3</sub> COO <input type="checkbox"/> +CF <sub>3</sub> COO	<input type="checkbox"/> H <sub>2</sub> O <input type="checkbox"/> H <sub>3</sub> PO <sub>4</sub>
<input type="text"/> <input type="button" value="X"/> <input type="button" value="+"/>	<input type="text"/> <input type="button" value="X"/> <input type="button" value="+"/>	<input type="text"/> <input type="button" value="X"/> <input type="button" value="+"/>

☐ Salt dominated positive ions (M+H may be weak or missing)

## Charge State Tab

Extraction	Ion Species	Charge State	Compound Filters	Mass Filters	Mass Defect	Results	Advanced
------------	-------------	--------------	------------------	--------------	-------------	---------	----------

Isotope grouping

Peak spacing tolerance:  m/z, plus  ppm

Isotope model: Common organic molecules

Charge state

☒ Limit assigned charge states to a maximum of:

☐ Treat ions with unassigned charge as singly-charged


## Compound Filters Tab

Extraction	Ion Species	Charge State	Compound Filters	Mass Filters	Mass Defect	Results	Advanced
<p>Height</p> <p><input checked="" type="checkbox"/> Relative height <math>\geq</math> 2.500 %</p> <p><input checked="" type="checkbox"/> Absolute height <math>\geq</math> 5000 counts</p> <p><input type="checkbox"/> Limit to the largest 100 compounds</p> <p>Compound location</p> <p><input type="checkbox"/> Restrict retention times to [ ] minutes</p> <p>Charge states</p> <p><input type="checkbox"/> Restrict charge states to [ ] Z</p>							

## Mass Filters Tab

Extraction	Ion Species	Charge State	Compound Filters	Mass Filters	Mass Defect	Results	Advanced
<p>Mass filters</p> <p><input checked="" type="checkbox"/> Filter mass list <math>\Delta</math> 5.000 ppm</p> <p>Include only these mass(es)</p> <p>Include only these mass(es)</p> <p>Exclude these mass(es)</p> <p><input checked="" type="radio"/> These masses:</p> <p>[ ]</p> <p>(type a comma-separated list of masses like "142.1012, 253.4003")</p> <p><input type="radio"/> Database</p> <p>D:\MassHunter\databases\default.csv [ ]</p>							

## Mass Defect Tab

Extraction	Ion Species	Charge State	Compound Filters	 Mass Filters	Mass Defect	Results	Advanced
------------	-------------	--------------	------------------	--	-------------	---------	----------

Mass defect filtering

☐ Filter results on mass defects

Expected mass defect

Constant

0.0000 Da + ( 0.0000 per 100.00 Da )


Calculate from formula

Mass defect tolerance

Constant (symmetric)

+/- 0.0100 Da

## Results Tab

Extraction	Ion Species	Charge State	Compound Filters	 Mass Filters	Mass Defect	Results	Advanced
------------	-------------	--------------	------------------	--	-------------	---------	----------

Previous results

☒ Delete previous compounds

New results

☒ Highlight first compound

☐ Highlight all compounds

Chromatograms and spectra

☐ Extract MFE spectrum ☐ Extract ECC

☐ Extract raw spectrum ☐ Extract EIC

☐ Prefer profile for raw spectrum, if available

☐ Clip extracted raw spectrum

Asymmetric (m/z) - 5.0000 + 10.0000

☐ Extract MS/MS Spectrum

☒ Separate MS/MS spectrum per CE

☐ Average MS/MS spectrum for all CEs

Precursor tolerance: +/- 20.00 ppm

☒ Deisotope MS/MS spectrum


Display limits

☐ Display only the largest 100 compounds


## Advanced Tab

Extraction	Ion Species	Charge State	Compound Filters	Mass Filters	Mass Defect	Results	Advanced
<p>Compound ion count threshold</p> <p> <input checked="" type="radio"/> Include all <input type="radio"/> Two or more ions <input type="radio"/> Only one ion </p>							
<p>Compounds with indeterminate neutral mass</p> <p> <input type="radio"/> Include <input checked="" type="radio"/> Exclude <input type="radio"/> Include only these </p>							

### 9.7.1.2 Method Explorer: Find Compounds by Formula-Find By Formula Options




## Formula Source Tab

 Formula Source   Formula Matching   Positive Ions   Negative Ions   Scoring   Results

Source of formulas to confirm

☐ These formulas:  
  
(type a comma-separated list of formulas, e.g., "C6H6, CH4")

☐ Compound exchange file (.CEF):

☒ Database   
D:\MassHunter\databases\default.csv

☐ Worklist


Values to match

☒ Mass

☐ Mass and retention time (retention time optional)

☐ Mass and retention time (retention time required)

## Formula Matching Tab

 Formula Source   Formula Matching   Positive Ions   Negative Ions   Scoring   Results

Match tolerance

Masses:                    +/- 5.00                    ppm

Retention times:                    +/- 0.50                    minutes

Expansion of values for chromatogram extraction

Possible m/z:   Symmetric (ppm)                     +/- 100.0

☒ Limit EIC extraction range

Expected retention time: +/- 1.50                    minutes






## Positive Ions Tab

Formula Source	Formula Matching	Positive Ions	Negative Ions	Scoring	Results
<div> <div> <p>Charge carriers</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> -electron</li> <li><input checked="" type="checkbox"/> +H</li> <li><input checked="" type="checkbox"/> +Na</li> <li><input type="checkbox"/> +K</li> <li><input type="checkbox"/> +NH4</li> </ul> </div> <div> <p>Neutral losses</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> H2O</li> </ul> </div> </div>					
<div> <div> <p>Charge states, if not known</p> <p>Charge state range 1-2</p> </div> <div> <p>Aggregates</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Dimers e.g., [2M+H]<sup>+</sup></li> <li><input type="checkbox"/> Trimers e.g., [3M+H]<sup>+</sup></li> </ul> </div> </div>					




## Negative Ions Tab

Formula Source	Formula Matching	Positive Ions	Negative Ions	Scoring	Results
<div> <div> <p>Charge carriers</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> +electron</li> <li><input checked="" type="checkbox"/> -H</li> <li><input type="checkbox"/> +Cl</li> <li><input type="checkbox"/> +Br</li> <li><input checked="" type="checkbox"/> +HCOO</li> <li><input type="checkbox"/> +CH3COO</li> <li><input type="checkbox"/> +CF3COO</li> </ul> </div> <div> <p>Neutral losses</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> H2O</li> </ul> </div> </div>					
<div> <div> <p>Charge states, if not known</p> <p>Charge state range 1-2</p> </div> <div> <p>Aggregates</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Dimers e.g., [2M-H]<sup>-</sup></li> <li><input type="checkbox"/> Trimers e.g., [3M-H]<sup>-</sup></li> </ul> </div> </div>					


## Scoring Tab

 Formula Source	Formula Matching	 Positive Ions	 Negative Ions	Scoring	Results
Contribution to overall score					
Mass score		<input type="text" value="100.00"/>			
Isotope abundance score		<input type="text" value="60.00"/>			
Isotope spacing score		<input type="text" value="50.00"/>			
Expected data variation					
MS mass:	<input type="text" value="2.0"/>	mDa	+	<input type="text" value="5.6"/>	ppm
MS isotope abundance:				<input type="text" value="7.5"/>	%
MS/MS mass:	<input type="text" value="5.0"/>	mDa	+	<input type="text" value="7.5"/>	ppm

## Results Tab

 Formula Source	Formula Matching	 Positive Ions	 Negative Ions	Scoring	Results
Previous results					
<input checked="" type="checkbox"/> Delete previous compounds					
New results					
<input checked="" type="radio"/> Highlight first compound					
<input type="radio"/> Highlight all compounds					
Matched results					
<input checked="" type="checkbox"/> Only generate compounds for matched formulas					
Chromatograms and spectra					
<input checked="" type="checkbox"/> Extract EIC					
<input checked="" type="checkbox"/> Extract raw spectrum					
<input type="checkbox"/> Extract cleaned spectrum					


### 9.7.1.3 Method Explorer: Identify Compounds-Search Database



### Search Criteria Tab


Search Criteria	Database	Peak Limits	Positive Ions	Negative Ions	Scoring	Search Mode	Search Results
<p>Values to match</p> <p> <input type="radio"/> Molecular formula  <input checked="" type="radio"/> Mass  <input type="radio"/> Mass and retention time (retention time optional)  <input type="radio"/> Mass and retention time (retention time required) </p> <p>Match tolerance</p> <p> Mass: <input type="text" value="5.00"/> <input type="text" value="ppm"/>  Retention time: <input type="text" value="0.10"/> <input type="text" value="minutes"/> </p>							

## Database Tab



One can choose a .csv, .mtl or a .cdb (PCDL file) here in dropdown menu.

## Peak Limits Tab



## Positive Ions Tab

Search Criteria	Database	Peak Limits	Positive Ions	Negative Ions	Scoring	Search Mode	Search Results
<div> <div> <p>Charge carriers</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> -electron</li> <li><input checked="" type="checkbox"/> +H</li> <li><input checked="" type="checkbox"/> +Na</li> <li><input type="checkbox"/> +K</li> <li><input type="checkbox"/> +NH4</li> </ul> </div> <div> <p>Neutral losses</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> H2O</li> </ul> </div> </div>							
<div> <div> <p>Charge states, if not known</p> <p>Charge state range: 1-2</p> </div> <div> <p>Aggregates</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Dimers e.g., [2M+H]<sup>+</sup></li> <li><input type="checkbox"/> Trimers e.g., [3M+H]<sup>+</sup></li> </ul> </div> </div>							



## Negative Ions Tab

Search Criteria	Database	Peak Limits	Positive Ions	Negative Ions	Scoring	Search Mode	Search Results
<div> <div> <p>Charge carriers</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> +electron</li> <li><input checked="" type="checkbox"/> -H</li> <li><input type="checkbox"/> +Cl</li> <li><input type="checkbox"/> +Br</li> <li><input checked="" type="checkbox"/> +HCOO</li> <li><input type="checkbox"/> +CH3COO</li> <li><input type="checkbox"/> +CF3COO</li> </ul> </div> <div> <p>Neutral losses</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> H2O</li> </ul> </div> </div>							
<div> <div> <p>Charge states, if not known</p> <p>Charge state range: 1-2</p> </div> <div> <p>Aggregates</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Dimers e.g., [2M-H]<sup>-</sup></li> <li><input type="checkbox"/> Trimers e.g., [3M-H]<sup>-</sup></li> </ul> </div> </div>							

## Scoring Tab

Search Criteria	Database	Peak Limits	Positive Ions	Negative Ions	Scoring	Search Mode	Search Results
<div> <p>Contribution to overall score</p> <p>Mass score: 100.00</p> <p>Isotope abundance score: 60.00</p> <p>Isotope spacing score: 50.00</p> </div>							
<div> <p>Expected data variation</p> <p>MS mass: 2.0 mDa + 5.6 ppm</p> <p>MS isotope abundance: 7.5 %</p> <p>MS/MS mass: 5.0 mDa + 7.5 ppm</p> </div>							

## Search Mode Tab

Search Criteria	Database	Peak Limits	 Positive Ions	 Negative Ions	Scoring	Search Mode	Search Results
-----------------	----------	-------------	---	---	---------	-------------	----------------

Ion search mode




Which database entries should be examined when searching masses from simple ions?

☐ Neutral entries


☒ Cation or anion entries

(This choice is not applicable to CSV databases.)

## Search Results Tab


Search Criteria	Database	Peak Limits	 Positive Ions	 Negative Ions	Scoring	Search Mode	 Search Results
-----------------	----------	-------------	---	---	---------	-------------	--

Search Results

☐ Limit to the best 

hits

### 9.7.1.4 Method Explorer: Identify Compounds-Generate Formulas



### Allowed Species Tab

Allowed Species Limits Charge State Scoring

Charge carrier to be assumed if not known

Positive ions: H Negative ions: H

MS ion electron state: allow both even and odd


Elements and limits

Element	Minimum	Maximum
C	3	60
H	0	120
O	0	30
N	0	10
S	0	5
Cl	0	10
F	0	50
P	0	5

+

X

## Limits Tab

 Allowed Species Limits Charge State Scoring

Limits on input masses

Maximum neutral mass for which formulas should be calculated: 1700.0000

Limits on results


☒ Minimum overall score 70.000

☐ Maximum MS mass error 7.5000 ppm

☐ Require DBE from 0.0 to 50.0

☐ Maximum number of hits 5

## Charge State Tab

 Allowed Species Limits Charge State Scoring

Isotope grouping

Peak spacing tolerance: 0.0025 m/z, plus 7.0 ppm

Isotope model: Common organic molecules

Charge state

☒ Limit assigned charge states to a maximum of: 2

☐ Treat ions with unassigned charge as singly-charged



## Scoring Tab

Allowed Species	Limits	Charge State	Scoring
<b>Contribution to overall score</b>			
Mass score			100.00
Isotope abundance score			60.00
Isotope spacing score			50.00
<b>Expected data variation</b>			
MS mass:	2.0	mDa +	5.6 ppm
MS isotope abundance:			7.5 %
MS/MS mass:	5.0	mDa +	7.5 ppm

## 9.8 Waste Management

All laboratory waste is disposed of as directed by the EPA-RTP Safety, Health, and Environmental Management Office.

### 9.8.1 Liquid Waste

Liquid waste generated during sample preparation is collected in a waste beaker during sample preparation and is emptied into a 4L empty solvent bottle that is designated and labeled for waste collection. Staff from Chemical Services are contacted and come to the lab to collect waste containers that are full.

### 9.8.2 Solid Waste

Solid waste generated during sample preparation is collected in beakers, plastic garbage bags, or burn boxes labeled for non-hazardous waste during sample preparation. All solid waste that is not collected in burn boxes during sample preparation is subsequently emptied into the burn box labeled for non-hazardous waste after sample preparation. Staff from Chemical Services is contacted to remove burn boxes that are full of solid laboratory waste.

## **10.0 Data and Records Management**

### **10.1 Laboratory Record Book**

Records of the preparation of samples, blanks, calibration standards, and QCs will be retained in a laboratory record book (LRB) that is kept by the individual conducting the analysis. This LRB will contain a record of all sample preparation activities and any other data that may be used to interpret results. All samples will be recorded in the LRB by a unique sample ID. The date of extraction and amount of internal standard/extraction solution made on each day of analysis will be recorded in the LRB. Data generated from analysis of samples will also be placed in the LRB. Generated instrumental data are regularly backed up on redundant storage space.

### **10.2 Location of LRB**

The LRB will be retained in the laboratory (or office area) where these operations are performed until the conclusion of the study and will be archived in a secure room for three years after completion of the study according to EPA records policy.

## **11.0 Quality Control and Quality Assurance**

Quality control samples are prepared with each batch of samples. Depending upon the assay being conducted this may be an SRM, spiked blank matrix or a sample that is simply run with every assay. This sample will be used to determine batch to batch consistency. As a measure of the accuracy of the analysis, the percent accuracy of the analytical value, relative to the pre-determined target value are calculated and reported along with the data as an indication of control when available. For NTA and suspect screening assays QA/QC samples are generally used for a different purpose. The purposes may include but are not limited to: batch to batch RT matching, as single point calibration of known compounds in the samples, for future anchoring with datasets not yet run. These QA/QC samples serves a variety of purposes, however it is not in the true sense the same as a QA/QC samples run in traditional analytical assays.

## **12.0 References**

1. EMAB-112.0- Extraction and Analysis of Perfluorinated Compounds (PFCs) from Surface Waters by High Performance Liquid Chromatography (HPLC)-Tandem Mass Spectrometry (MS/MS)
2. MDAB-23.0- Standard Operating Procedure (SOP) For Extraction And High Performance Liquid Chromatography/Mass Spectrometry Analysis Of Perfluorinated Acids And Sulfonates From House Dust

3. MDAB-077.0 - Extraction of Perfluorinated Compounds from Fish Liver Homogenate
4. PHCB-033-SOP-01 - Extraction and Analysis of SVOCs in Tire Crumb Rubber Samples
5. Agilent 6200 Series TOF LC/MS Techniques and Operation Course Number R1874A Volumes, I, II and III. Qual B.04.00 SP1; Quant B.04.00 SP3. Student Manual. July, 2011

Note: For copies of SOPs, please call the PI, Dr. Mark J Strynar (919-541-3706) or email at [Strynar.mark@epa.gov](mailto:Strynar.mark@epa.gov)