




TrichAnalytics Inc.



Method Validation for Trace Metals in Tissue by LA-ICP-MS

...and the reason this fish is jumping for joy!


Presented by: Jennie Christensen, PhD

TrichAnalytics Inc.

Outline

- Why LA-ICP-MS?
- Accuracy
- Precision
- Repeatability
- Sample Stability
- Detection Limits
- Internal Performance Testing
- Linearity/Range
- Robustness

 TrichAnalytics Inc.




Why LA-ICP-MS?

Current standard is bulk tissue sampling, ICP-MS, etc.

- Small sample volume = high detection limits
- High sample volume = lethality of individual, invasiveness, time
- >60,000 fish sacrificed since MMER of 2002 (Kambeitz et al. 2019)
- Wet to dry weight conversion error (high moisture tissues)
- Detail can be lost, some questions can't be answered with bulk analysis
- Mercury analyzed separately (cost)



TrichAnalytics Inc.



Why LA-ICP-MS?


Potential Benefits

- Small sample volume (<0.01 mm³)
- Do not consume sample, multiple analyses
- Invertebrates (more detail)
- Fish Tissue (non-lethal)
- Analyze samples dry – so no wet to dry weight conversion error
- Fast turnaround times
- No use of harsh chemicals
- Mercury included (cost benefit)



This sample has been
analyzed 25 times, so far!

(<1% volume consumed)



Trich Analytics Inc.

Why LA-ICP-MS?

Challenges to using LA-ICP-MS

- No standards (CRM) created for soft biological tissue for LA-ICP-MS
- Moisture in samples creates oxides on ablation
- Matrix-matching of standards and samples
- Potential analytical interference
- No similar laboratories for inter-lab testing
- No reference methods
- No lab currently accredited for LA-ICP-MS analysis of biological tissues (*to our knowledge)




 TrichAnalytics Inc.

We DID it! (or did we?) – enter, VALIDATION

TrichAnalytics Inc. – MET-002.02

Acknowledgements

- Incredible Trich team
- NAMC-SWG
- Teck Resources
- Allegiance Coal
- Minnow Environmental
- Nautilus Environmental
- Golder Associates



Accuracy

The measurement of how close the reportable (observed) result is to the true result.

$$\text{Accuracy (\%)} = \frac{\text{Reportable conc.}}{\text{True conc.}} \times 100$$

Method: DORM-4 (dogfish, muscle protein)

Test: NIST 1946 (wet, homogenized, fish muscle)

NIST 1566b (dry, powder, oyster)

NIST 2976 (dry, powder, mussel)



Trich Analytics Inc.

Precision

The measurement of variability of data from an experiment that has been repeated several times.

$$\text{Precision (\%)} = \frac{\text{St Dev of Reportable concs.}}{\text{Avg of Reportable concs}} \times 100$$

Method: DORM-4 (dogfish, muscle protein)

Test: NIST 1946 (wet, homogenized, fish muscle)

NIST 1566b (dry, powder, oyster)

NIST 2976 (dry, powder, mussel)



TrichAnalytics Inc.

Accuracy/Precision - Results

CRM	Observed	Known	Accuracy (%)	Precision (%)
DORM-4	3.35 ± 0.09	3.45 ± 0.40	97.1	2.6
NIST 1946	1.94 ± 0.16	1.72 ± 0.15	113	8.1
NIST 1566b	2.29 ± 0.08	2.06 ± 0.15	111	3.3
NIST 2976	2.59 ± 0.17	1.80 ± 0.15	144	6.7

N = 10 per CRM

Repeatability

Repeatability expresses the precision under the same operating conditions over a short interval of time – “intra-assay precision”

$$\text{Repeatability (\%)} = \frac{\text{St Dev all sample sets (ppm)}}{\text{Avg all sample sets (ppm)}} \times 100$$

Method: DORM-4 (dogfish, muscle protein)

Test: NIST 1946 (wet, homogenized, fish muscle)

NIST 1566b (dry, powder, oyster)

NIST 2976 (dry, powder, mussel)



Trich Analytics Inc.

Repeatability - Results

CRM	Observed	Actual	Repeatability (%)	Accuracy (%)
DORM-4	3.53 ± 0.12	3.45	3.5	102
NIST 1946	2.17 ± 0.10	1.715	8.4	127
NIST 1566b	2.33 ± 0.07	2.06	3.1	113
NIST 2976	2.58 ± 0.12	1.80	4.8	144

N = 30 per CRM

Sample Stability

The measurement of sample stability over time – change in concentration.

$$\text{Relative Percent Difference (\%)} = \frac{[\text{1}^{\text{st}} \text{ Reportable} - \text{2}^{\text{nd}} \text{ Reportable}]}{\text{1}^{\text{st}} \text{ Reportable}} \times 100$$

Method: DORM-4 (dogfish, muscle protein)

Test: 4 fish muscle samples (16 Aug 2019)


4 fish ovary samples (19 Jul 2019)

4 invertebrate samples (18-19 Jul 2019)



TrichAnalytics, Inc.

Sample Stability - Results



Excellent comparability
(slope = 1)

Little variability ($R^2 = 0.96$)

Large concentration range:
0.596 to 38.9 ppm

N=12

Detection Limits


Measurement of method detection limits to BC MOE (2016) Detection Quality Objectives for Aquatic Tissue (dry weight converted – 80% mc)

Element	BC MOE DQO	Element	BC MOE DQO
Lithium	-	Nickel	0.05
Sodium	10	Copper	0.05
Magnesium	10	Zinc	0.5
Aluminum	2	Arsenic	0.025
Phosphorus	25	Selenium	0.1
Potassium	50	Strontium	0.05
Calcium	10	Molybdenum	0.05
Vanadium	0.1	Cadmium	0.01
Chromium	0.05	Tin	0.1
Manganese	0.1	Mercury	0.01
Iron	5	Lead	0.02
Cobalt	0.02	Uranium	0.005



Trich Analytics Inc.

Detection Limits - Results



Internal Performance Testing

Comparison of reportable results using LA-ICP-MS to an independent CALA-certified laboratory (ICP-MS, etc.)

Linear regression – slope and r^2

Method: DORM-4 (dogfish, muscle protein)

Test: 12 benthic invertebrate samples


54 fish ovary samples

54 fish muscle samples



Trich Analytics, Inc.

Internal Performance Testing - Results



Linearity and Range

Measuring the relationship through calibration curve using standard linear regression analysis. And what range of concentrations can you accurately quantify?


Selenium – 4 point calibration using CRMs

160 point calibration using CRMs and samples



TrichAnalytics, Inc.

Linearity and Range - Results




Good linearity as per high R² value

But small range with just CRMs:
1.7 to 3.5 ppm

N=20

Linearity and Range - Results

Selenium



Excellent linearity as per high R^2 value and slope of 1

25X larger linear range:
0.6 to 47 ppm

N=160

Robustness

Consistency in results over time where machine variances occur (e.g., laser settings)

$$\text{Relative Percent Difference (\%)} = \frac{[\text{Modified} - \text{Unmodified}]}{\text{Average Result}} \times 100$$

Method: Standard laser settings (unmodified), DORM-4

Test: Laser Power x 2

Spot Size x 2

Laser Speed x 2



Robustness - Results

Parameter	Laser Speed	Laser Power	Spot Size
Unmodified	3.69 ± 0.08	3.83 ± 0.09	3.88 ± 0.12
Decrease	3.57 ± 0.07	3.53 ± 0.05	3.61 ± 0.17
Increase	3.60 ± 0.16	3.79 ± 0.16	3.86 ± 0.13

Parameter	Laser Speed	Laser Power	Spot Size
Unmodified	100%	100%	100%
Decrease	96.7%	91.8%	92.8%
Increase	97.5%	99%	99.5%

Summary

75% complete, 25 elements, submitting to CALA for accreditation
Preliminary Results: sample volume ~0.004 mm³

Validation Criteria	Selenium
Accuracy	Very Good (3/4 CRMs)
Precision	Excellent
Repeatability	Excellent
Sample Stability	Excellent
Detection Limits	Very Good (slightly above)
Internal Performance Testing	Excellent
Linearity and Range	Excellent
Robustness	Excellent



TrichAnalytics Inc.


Processing Speed

Manual (Qtegra and Excel) –
3 hours processing after
analysis

- Calibration slope
- Detection limits
- Concentration calculations
- Accuracy/Precision
- Duplicates


TrichAnalytics Inc.

Who you
callin' slow?



Processing Speed

In-House Data Automation/Programming:

3 hours to less than 1 min

- Calibration slope
- Detection limits
- Concentration calculations
- Accuracy/Precision
- Duplicates
- Diagnostic plots
- Processing log, including flags
- Output tables




TrichAnalytics Inc.

Told you I
was worth
hiring!



THANK YOU



TrichAnalytics Inc.