

Technical Memo

Shaw Environmental and Infrastructure, Inc.
Contract No. EP-C-08-034

In reply refer to: 12-JSC-04

Date: May 30, 2012

From: John Cox

To: Cindy Paul
Dennis Miller

Originator: Rick Wilkin

Thru: Sujith Kumar

Copies: Dominic Digiulio
Steve Vandegrift
Dennis Miller
John S. Cox

Subject: 8OA878DW

Task: 24274

Project: Geochemical Impacts

INTRODUCTION

Technical Memo 12-JSC-03 (TD8OA865CS) details the methodology used for the analysis of the 7 samples submitted under 8OA878DW. The results are to be considered as estimates as the C1-C4 Alcohol method is still in development and not all DQO's have been established.

METHOD AND MATERIALS

Column:

Restek RTX-624, 60 meter length, 0.53 mm ID with 3μ film thickness, catalog number 241580

OI Eclipse Purge-and-Trap with Infra-sparg heater:

Kit 4660 Sample heater 115v Item number 321732

Trap:

#7 Tenax only P&T Item number 227348

Standards:

Supelco C1-C4 Alcohols in water catalog number 21495418 Lot number LB90709 (methanol, ethanol, 2-propanol, tertiary butyl alcohol, n-propanol, 2-butanol, isobutanol, 1-butanol)

Restek C1-C4 Alcohols in water catalog number 566852 Lot number A086968

Working standards preparation:**C1-C4 Alcohols in water**

Supelco Lot# LB90709

units are ppb in 100 mls water

Stock A-10000 ppm in water (ampoule)

Stock B-100 ul of stock A in 900 uls water

ppm			ppb
Stock	amt of stock	amt water	Std Conc
1000 (B)	5 ul	100	50
1000 (B)	10 ul	100	100
1000 (B)	25 ul	100	250
10000(A)	5 ul	100	500
10000(A)	10 ul	100	1000

Eclipse Purge-and-Trap Method Modifications:

Only the parameters which were modified are listed. All other set points are as listed in RSKSOP-122 Revision 4 page 6.

Optimized P&T conditions for C1-C4 Alcohols in water		
Parameter	Standard	Optimized
Trap type	#10	#7
Sample temperature	Ambient	60 °C
Purge time	6 minutes	11 minutes
Desorb preheat temperature	off	On at 180°C
Trap temperature during purge	Ambient	40 °C
Water management fitting temperature during desorb	Ambient	40 °C

Instrument Method:

The instrument method is the same as that listed in RSKSOP-122 Revision 4, as seen below. The calibration has been changed to contain only C1-C4 alcohols and their isomers.

method: C:\CHEM32\2\METHODS\C1_C4_ALCOHOL.M\C1_C4_ALCOHOL.M
Modified on: 4/25/2012 at 9:49:13 AM

=====

6890 GC METHOD

=====

OVEN

Initial temp: 35 'C (On) Maximum temp: 260 'C
Initial time: 5.00 min Equilibration time: 0.20 min.
Ramps:
 # Rate Final temp Final time
 1 6.00 100 0.00
 2 20.00 200 8.00
 3 0.0(Off)
Post temp: 50 'C
Post time: 0.00 min
Run time: 28.83 min

FRONT INLET (SPLIT/SPLITLESS)

Mode: Split
Initial temp: 50 'C (Off)
Pressure: 0.00 psi (Off)
Total flow: 45.0 mL/min
Gas saver: Off
Gas type: Helium

BACK INLET (SPLIT/SPLITLESS)

Mode: Split
Initial temp: 150 'C (On)
Pressure: 6.74 psi (On)
Split ratio: 0.1:1
Split flow: 0.6 mL/min
Total flow: 10.1 mL/min
Gas saver: Off
Gas type: Helium

COLUMN 1

Capillary Column
Model Number: J&W 125-1364
DB-624
Max temperature: 260 'C
Nominal length: 60.0 m
Nominal diameter: 530.00 um
Nominal film thickness: 3.00 um
Mode: constant pressure
Pressure: 6.74 psi
Nominal initial flow: 5.1 mL/min
Average velocity: 33 cm/sec
Inlet: Back Inlet
Outlet: Front Detector
Outlet pressure: ambient

COLUMN 2

Capillary Column
Model Number: J&W 1211114
DB-624
Max temperature: 260 'C
Nominal length: 31.0 m
Nominal diameter: 320.00 um
Nominal film thickness: 1.80 um
Mode: (see column 1)
Pressure: 6.74 psi
Nominal initial flow: 1.3 mL/min
Average velocity: 23 cm/sec
Inlet: Back Inlet
Outlet: Back Detector
Outlet pressure: ambient

FRONT DETECTOR (CI4450)

Temperature: 250 'C (On)
Hydrogen flow: 40.0 mL/min (On)
Air flow: 375.0 mL/min (On)
Mode: Constant makeup flow
Makeup flow: 45.0 mL/min (On)
Makeup Gas Type: Helium
Flame: On
Electrometer: On
Lit offset: 2.0

BACK DETECTOR (PIDSIG)

Electrometer: On

SIGNAL 1

Data rate: 10 Hz

SIGNAL 2

Data rate: 20 Hz

method: C:\CHEM32\2\METHODS\C1_C4 ALCOHOL.M\C1_C4 ALCOHOL.M
Modified on: 4/25/2012 at 9:49:13 AM

Type: front detector
Save Data: On
Zero: 0.0 (Off)
Range: 0
Fast Peaks: Off
Attenuation: 0

Type: back detector
Save Data: Off
Zero: 0.0 (Off)
Range: 0
Fast Peaks: Off
Attenuation: 0

COLUMN COMP 1
Derive from front detector

COLUMN COMP 2
Derive from back detector

POST RUN
Post Time: 0.00 min

TIME TABLE

Time	Specifier	Parameter & Setpoint
------	-----------	----------------------

GC Injector

Front Injector:
No parameters specified

Back Injector:
No parameters specified

method: C:\CHEM32\2\METHODS\C1_C4_ALCOHOL.M\C1_C4_ALCOHOL.M
 Modified on: 4/25/2012 at 9:49:13 AM

Integration Events

Non signal specific Integration Events

Event	Value
Tangent Skim Mode	New Exponential
Tail Peak Skim Height Ratio	20.000
Front Peak Skim Height Ratio	20.000
Skim Valley Ratio	20.000
Baseline Correction	Advanced
Peak to Valley Ratio	500.000

Default Integration Event Table "Event"

Event	Value	Time
Initial Slope Sensitivity	1.000	Initial
Initial Peak Width	0.040	Initial
Initial Area Reject	1.000	Initial
Initial Height Reject	1.700	Initial
Initial Shoulders	OFF	Initial

Signal Specific Integration Event Table "Event_FID1A"

Event	Value	Time
Initial Slope Sensitivity	3.306	Initial
Initial Peak Width	0.103	Initial
Initial Area Reject	4.865	Initial
Initial Height Reject	0.395	Initial
Initial Shoulders	TAN	Initial

Apply Manual Integration Events: No

Calibration Table with Plots

```
=====
                        Calibration Table
=====
```

```
Calib. Data Modified :      4/24/2012 3:56:22 PM

Calculate           :      External Standard
Based on            :      Peak Area

Rel. Reference Window :      2.000 %
Abs. Reference Window :      0.000 min
Rel. Non-ref. Window :      2.000 %
Abs. Non-ref. Window :      0.000 min
Use Multiplier & Dilution Factor with ISTDs
Uncalibrated Peaks   :      not reported
Partial Calibration   :      Yes, identified peaks are recalibrated
Correct All Ret. Times:      No, only for identified peaks

Curve Type           :      Linear (some peaks differ, see below)
Origin               :      Included (some peaks differ, see below)
Weight               :      Equal (some peaks differ, see below)

Recalibration Settings:
Average Response      :      Floating Average New 0%
Average Retention Time:      Floating Average New 0%
```

```
Calibration Report Options :
  Printout of recalibrations within a sequence:
    Calibration Table after Recalibration
    Normal Report after Recalibration
  If the sequence is done with bracketing:
    Results of first cycle (ending previous bracket)
```

Signal 1: FID1 A,

RetTime [min]	Lvl Sig	Amount [ppb]	Area	Amt/Area	Ref Grp Name
6.399	1 1	50.00000	7.72840	6.46965	Methanol
	2	100.00000	14.37487	6.95659	
	3	250.00000	33.48655	7.46568	
	4	500.00000	65.95109	7.58138	
	5	1000.00000	127.33865	7.85308	
	6	2000.00000	253.71498	7.88286	
8.391	1 1	50.00000	60.17266	8.30942e-1	Ethanol
	2	100.00000	118.34559	8.44983e-1	
	3	250.00000	305.34064	8.18758e-1	
	4	500.00000	602.22833	8.30250e-1	
	5	1000.00000	1192.22107	8.38771e-1	
	6	2000.00000	2349.09814	8.51391e-1	
9.843	1 1	50.00000	173.56769	2.88072e-1	2-Propanol
	2	100.00000	328.54593	3.04371e-1	
	3	250.00000	851.23224	2.93692e-1	
	4	500.00000	1712.65564	2.91944e-1	
	5	1000.00000	3338.23096	2.99560e-1	
	6	2000.00000	6718.08301	2.97704e-1	
10.938	1 1	50.00000	436.51917	1.14543e-1	Tert Butyl Alcohol
	2	100.00000	845.62994	1.18255e-1	

RetTime [min]	Lvl Sig	Amount [ppb]	Area	Amt/Area	Ref Grp Name
		3 250.00000	2181.03369	1.14625e-1	
		4 500.00000	4295.52930	1.16400e-1	
		5 1000.00000	8465.59375	1.18125e-1	
		6 2000.00000	1.69248e4	1.18170e-1	
12.677	1	1 50.00000	159.04778	3.14371e-1	n-propanol
		2 100.00000	301.40417	3.31780e-1	
		3 250.00000	786.15668	3.18003e-1	
		4 500.00000	1584.15601	3.15625e-1	
		5 1000.00000	3054.05835	3.27433e-1	
		6 2000.00000	6170.49902	3.24123e-1	
14.365	1	1 50.00000	317.38031	1.57540e-1	2-Butanol
		2 100.00000	608.55652	1.64323e-1	
		3 250.00000	1569.97827	1.59238e-1	
		4 500.00000	3145.65894	1.58949e-1	
		5 1000.00000	6179.34570	1.61829e-1	
		6 2000.00000	1.24065e4	1.61205e-1	
15.762	1	1 50.00000	391.51703	1.27708e-1	Isobutanol
		2 100.00000	751.86346	1.33003e-1	
		3 250.00000	1931.40845	1.29439e-1	
		4 500.00000	3878.64111	1.28911e-1	
		5 1000.00000	7664.20996	1.30477e-1	
		6 2000.00000	1.53577e4	1.30228e-1	
17.198	1	1 50.00000	253.30080	1.97394e-1	1-Butanol
		2 100.00000	475.53476	2.10290e-1	
		3 250.00000	1242.83899	2.01152e-1	
		4 500.00000	2474.64551	2.02049e-1	
		5 1000.00000	4787.57373	2.08874e-1	
		6 2000.00000	9630.36133	2.07677e-1	

More compound-specific settings:

Compound: Methanol
Curve Type : Linear
Origin : Included
Calibration Level Weights:/
Level 1 : 1
Level 2 : 1
Level 3 : 1
Level 4 : 1
Level 5 : 1
Level 6 : 1

Compound: Ethanol
Curve Type : Linear
Origin : Included
Calibration Level Weights:/
Level 1 : 1
Level 2 : 1
Level 3 : 1
Level 4 : 1
Level 5 : 1
Level 6 : 1

Compound: 2-Propanol
Curve Type : Linear
Origin : Included

Grendel 5/23/2012 12:39:01 PM jsc

Page 2 of 5

Method C:\CHEM32\2\METHODS\C1_C4_ALCOHOL.MAC1_C4_ALCOHOL.M

Compound: Tert Butyl Alcohol
Curve Type : Linear
Origin : Included

Compound: n-propanol
Curve Type : Linear
Origin : Included





Compound: 2-Butanol
Curve Type : Linear
Origin : Included
Calibration Level Weights:/
Level 1 : 1
Level 2 : 1
Level 3 : 1
Level 4 : 1
Level 5 : 1
Level 6 : 1





Compound: Isobutanol
Curve Type : Linear
Origin : Included

Compound: 1-Butanol
Curve Type : Linear
Origin : Included

=====

Calibration Curves





Sample Analysis:

The test samples were analyzed using the method outlined above. Several blanks were analyzed at the beginning of the sample queue followed by a Continuing Calibration Check standard (CCC) at 100 ppb and a Second Source Standard (SS) at 500 ppb.

The samples were analyzed sequentially following the Blank, CCC, and SS. The results of the analysis are summarized on table 1.

Table 1: Sample analysis results-analysis date 15-May-2012

80A878DW

BLACK= FIRST RUN

units

WITH FULL 40 ML

=ppb

VOA

SS 6359

RED=SECOND RUN WITH 5ML HEADSPACE USING
SAME VOA VIALS AS FIRST RUN

SAMPLE	METHANOL	ETHANOL	2- PROPANOL	TERT BUTYL ALCOHOL	N- PROPANOL	2- BUTANOL	ISOBUTANOL	1- BUTANOL
PGDW20-0412	ND	ND	ND	ND	ND	ND	ND	ND
PGDW20-0412	ND	ND	ND	ND	ND	ND	ND	ND
EPAMW02-0412	592	28.4	51.7	5395	11.8	15.6	10.7	32.2
EPAMW02-0412	337	28.7	42.7	5448	10.6	10.4	10.9	29.6
PGDW30-0412	ND	ND	ND	ND	ND	ND	ND	ND
PGDW30-0412	ND	ND	ND	BQL-1.2	ND	ND	ND	ND
PGDW05-0412	ND	ND	ND	ND	ND	ND	ND	ND
PGDW05-0412	ND	ND	ND	ND	ND	ND	ND	ND
EPAMW01-0412	863	ND	ND	197	ND	ND	BQL-5.8	ND
EPAMW01-0412	894	ND	ND	54.4	ND	ND	BQL-2.1	ND
EPAMW01-0412D	831	ND	BQL-3.9	198	BQL-1.7	ND	BQL-5.5	BQL-1.9

EPAMW01-0412D	837	BQL-13.7	BQL-1.1	60	ND	ND	BQL-1.8	BQL-2.6
FIELDCLK4	ND	ND	ND	ND	ND	ND	ND	ND
FIELDCLK4	ND	ND	ND	ND	ND	ND	ND	ND

RESULTS AND DISCUSSION

The VOA vials containing the CCC and the SS were prepared as duplicates to those utilized in the ICAL (24-April-2012). The CCC was within 1% for all alcohols except methanol and ethanol with the ethanol value at -7% and the methanol at +200%. The SS was within 5% of the true value. The ending CCC was done at both 50ppb and 250ppb with recoveries in the 70-130% range. Subsequent reanalysis of freshly prepared CCC's demonstrated the validity of the calibration curve with recovery error less than 10% for all alcohols. This data demonstrates there are problems with refrigerated storage of prepared standards. In this case, the vials were stored at 4°C for about 3 weeks.

There are two sets of data in Table 1; the data in black is the first run of the VOA vials without any headspace. The second run of the vials (in red) was done the second day after the samples had set on the autosampler, at ambient temperature, overnight with 5 mls of headspace. There is one field duplicate in the sample set, ss#6359-5 and 6359-6. The RPD for the methanol value is 0.72% and the same value for TBA is 0.51%.

Literature Cited:

“Low-level Detection of Ethanol, 1,4-Dioxane, and Other Oxygenates Using the Eclipse Purge-and-Trap Sample Concentrator”. OI Analytical Application Note 25370206 12 pages. March 2006.

RSKSOP-122/Revision 4. Analysis of Volatile Aromatic Hydrocarbons by Purge and Trap Gas Chromatography. John Cox. March 2007.

RSKSOP-248/Revision 1, Standard Operating Procedure Determination of Method Detection Limits. Steve Vandegrift, May 2004.

EPA method 8015C Revision 3. Nonhalogenated Organics by Gas Chromatography. February 2007.