

TEST METHOD

RESTRICTED

SUBJECT DETERMINATION OF AREA PERCENT PURITY OF EASTMAN™ PRODUCTS, CHDM-D, VERITYL AB-1000, CHDM-D90 GLYCOLS AND HIGH-BOILING CHDM DISTILLATES BY GAS CHROMATOGRAPHY	SEQUENCE NUMBER
	EFFECTIVE DATE September 22, 2011
	PAGE 1 OF 6 PAGES
	ACCESSION NUMBER
	ECCN EAR99
METHOD NO. EMN-A-AC-G-GC-500-7-SPEC	SUPERSEDES ED. 6 DATED 10/28/2009

EXCEPT FOR DESIGNATED CONTROLLED COPIES, A PRINTED HARD COPY IS VALID ONLY 24 HOURS PAST 1/11/14 4:37 PM.

1. SCOPE AND SIGNIFICANCE

This method describes a gas chromatographic procedure for determination of the area percent purity of 1,4-cyclohexanedimethanol (CHDM) and reporting the percent of the major impurities present in CHDM-D, Verityl AB-1000, and CHDM-D90 glycols. This method is applicable to all CHDM samples from the refining columns. The range of this method on any of the impurities is from 0.05 to 50 percent, and for the CHDM is 0.1 to 99.9 percent. The analyst time is ten minutes, and the instrument time is twenty-five minutes.

2. SUMMARY OF METHOD

The sample is dissolved in dimethylformamide (DMF) and chromatographed on a (50%-phenyl)-methylpolysiloxane capillary column using a flame ionization detector. The concentrations of sample components are calculated from the integrated chromatogram using area normalization.

3. DEFINITIONS

The terms employed in this method are commonly used in normal laboratory practice and require no special comment.

4. PRECISION AND BIAS

4.1 Precision

The relative standard deviations are based on thirty injections of one prepared standard. Relative standard deviations ranged from 0.31% to 1.13% for individual components.

4.2 Bias

No calibration standards are used for area percent determinations; therefore, the accuracy has not been determined.

For reasons of safety and accuracy, the person performing this procedure must be thoroughly trained and under the supervision of a professional person who is knowledgeable in the relevant science. Equipment and materials described should be used in accordance with safety precautions recommended by their manufacturers.

Neither Eastman Chemical Company nor its marketing affiliates shall be responsible for the use of this information, or of any product, method, or apparatus mentioned, and you must make your own determination of its suitability and completeness for your own use, for the protection of the environment, and for the health and safety of your employees and purchasers of your products. No warranty is made of the merchantability or fitness of any product; and nothing herein waives any of the seller's conditions of sale.

EASTMAN

**Eastman Chemical Company
Kingsport, Tennessee 37662**

5. SAFETY PRECAUTIONS

- 5.1 Nitrile gloves must be worn when handling DMF. Work should be performed in a well-ventilated area to minimize vapor inhalation. If exposure to the vapors cannot be avoided (for example, if the sample must be heated without a cover or volumes greater than 20 mL must be transferred), the work must be performed in a hood.
- 5.2 Otherwise, normal safety precautions and safe handling practices should be observed.
- 5.3 Refer to the PPE Matrix posted in the lab for required personal protective equipment.

6. SOURCES OF ERROR

The method is an area percent determination where each component is given the same detector response factor of 1.0000, and the total data is normalized to 100%. If the assumption of equal detector response factors for every component is incorrect, the area percent cannot be taken to represent the weight percent of each component in the solution. Also, the method will not account for non-eluting components.

7. APPARATUS

- 7.1 Gas chromatograph with split injector and FID
- 7.2 Autosampler
- 7.3 Chromatography data system
- 7.4 Vials, autosampler with septum caps
- 7.5 Capillary column, fused-silica, (50%-phenyl)-methylpolysiloxane, 30-m x 0.32-mm id with film thickness = 0.25 μ m
- 7.6 Disposable transfer pipets
- 7.7 Vials, 8-dram
- 7.8 Karl Fischer titrator

8. REAGENTS AND MATERIALS

Dimethylformamide (DMF) (R-3, S-3, F-2, C-0)

9. CALIBRATION AND STANDARDIZATION

This method uses an area percent calculation because of the difficulty in obtaining pure materials. If possible, a sample of crude CHDM containing high levels of impurities should be obtained and used to verify the retention times.

10. SAMPLE PREPARATION

- 10.1 Ten grams of the colorless sample are required for a single analysis.

- 10.2 If the sample as received has two layers, shake vigorously before taking a sample. If the sample has solidified, warm it up in an oven at approximately 100°C until it is a liquid.

11. PROCEDURE

- 11.1 If the sample label asks for percent water by Karl Fischer titration, the GC results are corrected for water automatically in LIMS.
- 11.2 Transfer approximately 0.2 grams of sample into an 8-dram vial.
- 11.3 Add 5 mL of dimethylformamide (DMF) to the vial, cap and shake well to dissolve.
- 11.4 Transfer approximately 1 mL of the sample prepared in Step 11.3 to an autosampler vial and cap the vial.
- 11.5 Record and integrate the chromatogram using the conditions specified in Section 7 and Appendix A1.

12. CALCULATIONS

All calculations are performed by the chromatography data system, or the area percents of each component can be calculated by hand as follows:

$$\text{Area \% Component} = \frac{(\text{Area of Component}) \times 100}{\text{Sum of Areas of All Detected Components Excluding DMF Solvent}}$$

$$\text{Area \% Component (corrected for water)} = \frac{(\text{Area of Component}) \times (100 - \% \text{ Water})}{\text{Sum of Areas of All Detected Components Excluding DMF Solvent}}$$

$$\text{ppm} = \text{Area\%} \times 10,000$$

13. REPORT

Report the percents of the identified components to one decimal place.

14. BIBLIOGRAPHY

None

15. APPENDIX

A1. Instrument Conditions

INSTRUMENT CONDITIONS

(Page 1 of 3)

Oven Parameters

Oven: On
 Oven Equilib Time: 0.50 min
 Oven Max. Temperature: 320 C
 Oven Cyro: Off
 Oven Cyro Blast: Off
 Oven Ambient: 20 C
 Timeout Detection: Off
 Timeout Detector: 10.00 min
 Fault Detection: Off

Oven Temperature Program

Level	Rate	Next Temp	Hold Time
Initial		115.00 C	0.00 min
1	5C/min	130.00 C	0.00 min
2	100 C/min	200.00 C	0.00 min
3	120 C/min	280.00 C	3.00 min

Zone Temperatures

Front Inlet Setpoint: 280 C
 Back Inlet Setpoint: 280 C
 Front Detector Setpoint: 280 C
 Back Detector Setpoint: 280 C

Inlet Temperature Setup

Inlet: Front

Oven Track: Off
 Zone Temperature: 280 C
 Inlet Temperature Program: Off

Inlet: Back

Oven Track: Off
 Zone Temperature: 280 C
 Inlet Temperature Program: Off

Inlet Pressure Setup

Inlet: Front

General Inlet Settings
 S/SL, COC, PP Inlet

Split Mode: Split
 Split Flow: 269 mL/min
 Split Ratio: 75.00 : 1
 Gas Saver: 20.00 mL/min @ 2.00 min

INSTRUMENT CONDITIONS
(Page 2 of 3)

Inlet: Back

General Inlet Settings

S/SL, COC, PP Inlet

Split Mode:	Split
Split Flow:	269 mL/min
Split Ratio:	75.00 : 1
Gas Saver:	20.00 mL/min @ 2.00 min

Column Setup

Column 1

Capillary Column

Length:	30.00 m
Inside Diameter:	0.25 mm
Film Thickness:	0.25 um
Inlet:	Front
Detector:	Front
Gas:	Hydrogen

Carrier Flow

Column Mode:	Constant Flow
Flow:	3.50 mL/min

Column 2

Capillary Column

Length:	30.00 m
Inside Diameter:	0.25 mm
Film Thickness:	0.25 um
Inlet:	Back
Detector:	Back
Gas:	Hydrogen

Carrier Flow

Column Mode:	Constant Flow
Flow:	3.50 mL/min

Signals Setup

Signal 1 Parameters

Detector:	Front Det
Range:	0
Attenuation:	0
Auto Zero:	Off

Signal 2 Parameters

Detector:	Back Det
Range:	0
Attenuation:	0
Auto Zero:	Off

INSTRUMENT CONDITIONS

(Page 3 of 3)

Run Time Control Setup

Run Time Control Program

#	Time (min)	Event	Parameter	Value
Front FID Setup				
		Heater:	280 C	
		H2 Flow:	30 mL/min	
		Air Flow:	350 mL/min	
		Makeup Flow:		
		Makeup Flow Type:	Helium	
		Constant Column:	30 mL/min	
		Flame:	On	
		Electrometer:	On	

Back FID Setup

	Heater:	280 C
	H2 Flow:	30 mL/min
	Air Flow:	350 mL/min
	Makeup Flow:	
	Makeup Flow Type:	Helium
	Constant Column:	30 mL/min
	Flame:	On
	Electrometer:	On

Autosampler Instrument Setup

Front Injector Program

Sample Washes:	4
Sample Pumps:	4
Viscosity Delay:	0
Solvent A Post-washes:	4
Solvent B Post-washes:	4

Rear Injector Program

Sample Washes:	4
Sample Pumps:	4
Viscosity Delay:	0
Solvent A Post-washes:	4
Solvent B Post-washes:	4

Trigger Type: External Trigger