Data Validation SOP HW-2, Rev. 11 Metals

Evaluation of Metals Data for the Contract Laboratory Program (CLP) based on

	SOW.	3/90
_		

(SOP Revision XI)

PREPARED	BY:		DATE:
		Hanif Sheikh, Quality Assurance Chemist Toxic and Hazardous Waste Section	
APPROVED	BY:	Kevin Kubik, Chief Toxic and Hazardous Waste Section	DATE:
APPROVED		Robert Runyon, Chief Monitoring Management Branch	DATE:

STANDARD OPERATING PROCEDURE

Page 1 of 34

Title: Evaluation of Metals Data for the Contract Laboratory Program

Date: Jan. 1992 Number: HW-2

Revision: 11

1.0 Scope

- This procedure is applicable to inorganic data obtained from contractor 1.1 laboratories working for Hazardous Waste Site Contract Laboratory Program (CLP).
- The data validation is based upon analytical and quality assurance 1.2 requirements specified in Statement of Work (SOW) 3/90 .
- 2.0 Responsibilities - Data reviewers will complete the following tasks as assign the Data Review Coordinator:
 - 2.1. For a total review:
 - 2.1.1 Data Assessment "Total Review-Inorganics" Checklist Appendix (A.1). The reviewer must answer every question on the checklist.
 - 2.1.2 Data Assessment Data Assessment Narrative (Appendix A.2) The answer on the checklist must match the action in the narrative (appendix A.2) and on Form I's. Do not use pencil to write the narrative.
 - 2.1.3 Contract Non-Compliance SMO Report (Appendix A.3) This report is to be completed only when a serious contract violation is encountered, or upon the request of the Data Validation Task Monitor, or Tech Project Officer (TPO). Forward 5 copies: one each for internal files, appropriate Regional TPO, Sample Management Office (SMO) and last two address Mailing List for Data Reviewers (Appendix A.4). In other cases, all contract violations should be appended to the end of the Data Assessment Narrative (Se A.2.2).

2.1.4 CLP Data Assessment Summary Forms

2.1.4.1 Appendix A.5

Fill in the total number of analytes analyzed by different analyses and the number of analytes rejected or flagged as estimated due to corresponding quality control criteria. Place an "X" in boxes where analyses were not performed, or criteria do not apply.

2.1.4.2 Appendix A.6

Data reviewer is also required to fill out Inorganic Regional Data Assessment form (Appendix A.7) provided by EPA Headquarters. Codes listed on the form will be used to describe the Data Assessment Summary.

STANDARD OPERATING PROCEDURE

Page 2 of 34

Title: Evaluation of Metals Data for the Contract Laboratory Program Date: Jan. 1992 Number: HW-2 Revision: 11

- 2.1.5 **Data Review Log:** It is recommended that each data reviewer should maintain a the reviews completed to include: a. date of start of case
 - b. date of completion of case review
 - c. site
 - d. case number
 - e. contract laboratory
 - f. number of samples
 - g. matrix
 - h. hours worked
 - i. reviewer's initials
- 2.1.6 Telephone Record Log the data reviewer should enter the bare facts of inquiry, before initiating any phone conversation with CLP laboratory. After the case review has been completed, mail white copy of Telephone Record Log to the laboratory and pink copy to SMO. File yellow copy in the Telephone Record Log folder, and attach a xerox copy of the Telephone Record Log to the completed Data Assessment Narrative (Appendix A.2).

2.1.7 Forwarded Paperwork

- 2.1.7.1 Upon completion of review, the following are to be forwarded to the Regional Sample Control Center (RSCC) located in the Surveillance and Monitoring Branc
 - a. data package
 - b. completed data assessment checklist (Appendix A.1, original)
 - c. SMO Contract Compliance Screening (CCS)
 - d. Record of Communication (copy)
 - e. CLP Reanalysis Request/Approval Record (original + 3 copies)
 - f. Appendix A.6 (original).
- 2.1.7.2 Forward 2 copies of completed Data Assessment Narrative (Appendix A.2) along with <u>2</u> copies of the Inorganic Data Assessment Form (Appendix A.6) and Telephone Record Log , if any,: one each for appropriate Regional TPO, and the other one to EPA EMSL office in Las Vegas. The addresses of TPOs and office in Las Vegas are given in Appendix A-4.
- 2.1.8 **Filed Paperwork** Upon completion of review, the following are to be filed within MMB files:
 - a. Two copies of completed Data Assessment Narrative (Appendix A.2) each carr Appendix A.6.
 - b. Telephone Record Log (copy)
 - c. SMO Report (copy Appendix A-3)
 - d. CLP Reanalysis Request/Approval Record (copy)

STANDARD OPERATING PROCEDURE

Page 3 of 34

Title: Evaluation of Metals Data for the Contract Laboratory Program

Date: Jan. 1992 Number: HW-2 Revision: 11

3.0 <u>Data Completeness</u>

Each data package is checked by a Regional Sample Control Coordinator (RSSC) completeness. A data package is assumed to be complete when all the delivera required under the contract are present. If a data package is incomplete, the would call the laboratory for missing document(s). If the laboratory does not within a week, SMO and MMB coordinator of Region II will be notified

4.0 Rejection of Data - All values determined to be unacceptable on the Inorganic Analysis Data Sheet (Form I) must be lined over with a red pencil. As soon review criteria causes data to be rejected, that data can be eliminated any further review or consideration.

Acceptance Criteria - In order that reviews be consistent among reviewers, acceptance criteria as stated in Appendix A.1 (pages 4-25) should be used. Additional guidance can be found in the National Inorganic Functional Guidel October 1, 1989. 6.0 SMO Contract Compliance Screening (CCS) - This is intended to aid reviewer in locating any problems, both corrected and uncorrected. However, the validat should be carried out even if CCS is not present. Resubmittals received from laboratory in response to CCS must be used by the reviewer. Request for Reanalysis - Data reviewers must note all items of contract non-compliance within Data Assessment Narrative. If holding times and sample times have not been exceeded, TPO may request reanalysis if items of non-com are critical to data assessment. Requests are to be made on "CLP Re-Analysis Request/Approval Record". 8.0 Record of Communication - Provided by the Regional Sample Control Center (RSC indicate which data packages have been received and are ready to be reviewed. 9.0 Rounding off numbers - The data reviewer will follow the standard practice. STANDARD OPERATING PROCEDURE Page 4 of 34 Title: Evaluation of Metals Data for the Date: Jan. 1992 Contract Laboratory Program Number: HW-2 Appendix A.1: Data Assessment - Contract Revision: 11 Compliance (Total Review) YES NO A.1.1 Contract Compliance Screening Report (CCS) - Present? [] ACTION: If no, contact RSCC. A.1.2 Record of Communication (from RSCC) - Present? [] ACTION: If no, request from RSCC.

A.1.6	Form I to IX	Yes	No
Title:	Evaluation of Metals Data for the Contract Laboratory Program Appendix A.1: Data Assessment - Contract Compliance (Total Review)	Date: Ja Number: Revision:	HW-2
	STANDARD OPERATING PROCEDURE	Page 5	of 34
	ACTION: If no for any of the above, contact RSCC for clarification.		
	(b) Form I's?	[]	
	(a) Traffic Report Sheet?	[]	
	Do sample numbers on cover page agree with sample numbers on:		
	Do numbers of samples correspond to numbers on Record of Communication?	[]	
	ACTION: If no, prepare Telephone Record Log, and contact laboratory.		
	Is cover page properly filled in and signed by the lab manager or the manager's designee?	[]	
A.1.5	Cover Page - Present?	[]	-
	ACTION: If no, request from Regional Sample Control Center (RSCC).		
	Legible?	[]	
A.1.4	Sample Traffic Report - Present?	[]	
	ACTION: If no, contact RSCC for trip report.		
A.1.3	Trip Report - Present and complete?	[]	

A.1.6.1	Are all the Form I through Form IX labeled with:		
	Laboratory	name? []	
	Case/SAS nu	mber? []	
	EPA sample	No.? []	
	SDG	No.? []	
	Contract	No.? []	
	Correct us	nits? []	
	Mat	trix? []	
	ACTION: If no for any of the above, note under Contract Problem/Non-Compliance section of the "Data Assessment Narrative".		
A.1.6.2	Do any computation/transcription errors exceed 10 reported values on Forms I-IX for:)% of	
	(NOTE: Check all forms against raw data.)		
	(a) all analytes analyzed by ICP?	[]	
	(b) all analytes analyzed by GFAA?	[]	
	(c) all analytes analyzed by AA Fla	ame? []	
	(d) Mercury?	[]	
	(e) Cyanide?	[]	
	ACTION: If yes, prepare Telephone Log, contact laboratory for corrected data and correct errors with red pencil and initial	al.	

Appendix A.1: Data Assessment - Contract Revision: 11 Compliance (Total Review) YES NO A.1.7 Raw Data A.1.7.1 Digestion Log* for flame AA/ICP (Form XIII) present? Digestion Log for furnace AA Form XIII present? [] Distillation Log for mercury Form XIII present?] Distillation Log for cyanides Form XIII present? Are pH values (pH<2 for all metals, pH>12 for cyanide) present? *Weights, dilutions and volumes used to obtain values. Percent solids calculation present for soils/sediments? [] Are preparation dates present on sample preparation logs/bench sheets? [] A.1.7.2 Measurement read out record present? ICP [] Flame AA Furnace AA [] Mercury [] Cvanides A.1.7.3 Are all raw data to support all sample analyses and QC operations present? [] Legible? [] Properly Labeled? [] ACTION: If no for any of the above questions

in sections A.1.7.1 through A.1.7.3, write Telephone Record Log and contact

Date: Jan. 1992

HW-2

Number:

Title: Evaluation of Metals Data for the

Contract Laboratory Program

laboratory for resubmittals.

	STANDARD OPERATING PROCEDURE	Page	7 of 34
Title:	Evaluation of Metals for the Contract Laboratory Program Appendix A.1: Data Assessment - Contract	Number	Jan. 1992 : HW-2 on: 11
	Compliance (Total Review)		
A.1.8	Holding Times - (aqueous and soil samples)	YES	NO
	(Examine sample traffic reports and digestion/distilla-	tion logs.	.)
	Mercury analysis (28 days) exceeded?		[] _
	Cyanide distillation (14 days) exceeded?		[]
	Other Metals analysis (6 months) exceeded?		[]
	NOTE: Prepare a list of all samples and analytes for which holding times have been exceeded. Specthe number of days from date of collection to of preparation (from raw data). Attach to chemical days from the collection of preparation (from raw data).	cify o the date	<u> </u>
	ACTION: If yes, reject (red-line) values less than Instrument Detection Limit (IDL) and flag as estimated (J) the values above IDL even though sample(s) was preserved properly.		
A.1.8.2	Is pH of aqueous samples for: Metals Analysis >2?		[]
	Cyanides Analysis <12?		[]
	Action: If yes, flag the associated metals and cyanic data as estimated.	des	
A.1.9	Form I (Final Data)		
A.1.9.1	Are all Form I's present and complete?	[]	

laboratory for submittal. A.1.9.2 Are correct units (ug/l for waters and mg/kg for soils) indicated on Form I's? Are soil sample results for each parameter corrected for percent solids? Are all "less than IDL" values properly coded with "U"? [STANDARD OPERATING PROCEDURE Page of 34 Title: Evaluation of Metals Data for the Date: Jan. 1992 Contract Laboratory Program Number: HW-2Appendix A.1: Data Assessment - Contract Revision: 11 Compliance (Total Review) YES NO Are the correct concentration qualifiers used with final data? If no for any of the above, prepare Telephone ACTION: Record Log, and contact laboratory for corrected data. A.1.9.3 Are EPA sample # s and corresponding laboratory sample ID # s the same as on the Cover Page, Form I's and in the raw data? [] Was a brief physical description of samples given on Form I's? Was the dilution of any sample diluted beyond the requirements of the contract noted on Form I or Form XIV? ACTION: If no for any of the above, note under Contract-Problem/Non-Compliance of the "Data Assessment Narrative". A.1.10 Calibration

ACTION: If no, prepare telephone record log and contact

A.1.10.1 Is record of at least 2 point calibration

present for ICP analysis?	[]	
Is record of 5 point calibration present for Hg analysis?	[]	
Is record of 4 point calibration present for:		
Flame AA?	[]	
Furnace AA?	[]	
Cyanides?	[]	
Is one calibration standard at the CRDL level for all AA (except Hg) and cyanides analyses?	[]	
ACTION: If no for any of the above, write in the Contract Problem/Non-Compliance section of the "Data Assessment Narrative".		
STANDARD OPERATING PROCEDURE	Page	9 of 34
Title: Evaluation of Metals Data for the Contract Laboratory Program Appendix A.1: Data Assessment - Contract Compliance (Total Review)	Date: Number Revisi	
A.1.10.2 Is correlation coefficient less than 0.995 for:	YES	NO
Mercury Analysis?		[]
Cyanide Analysis?		[]
Atomic Absorption Analysis?		[]
ACTION: If yes, flag the associated data as estimated.		
NOTE: The data validator shall calculate the correlation coefficient using concentrations of the standards and the corresponding instrument response (e.g. absorbance, peak area, peak height, etc.).	3	
A.1.10.3 In the instance where less than 4 standards are measured in absorbance (or peak area, peak height mode, are the remaining standards analyzed in	,etc.)	

	ACTION: If no, flag the associated data as estimated if standards are not within ±10% of true values Do not flag the data as estimated in linear ran indicated by good recovery of standard(s).	
A.1.11	Form II A (Initial and Continuing Calibration Verification	.on) -
A.1.11.1	Present and complete for every metal and cyanide?	[]
	Present and complete for AA and ICP when both are used for the same analyte?	[]
	ACTION: If no for any of the above, prepare Telephone Record Log and contact laboratory.	
A.1.11.2	Circle on each Form IIA all percent recoveries that are outside the contract windows. Are all calibration standards (initial and continuing) within control limits:	
	Metals- 90-110%R?	[]
	Hg - 80-120%R?	[]
	Cyanides- 85-115%R?	[]
	STANDARD OPERATING PROCEDURE	Page 10 of 34
C P	Evaluation of Metals Data for the Contract Laboratory Program Appendix A.l: Data Assessment - Contract Compliance (Total Review)	Date: Jan. 1992 Number: HW-2 Revision: 11
ACT	ION: Flag as estimated (J) all positive data (not flagged with a "U") analyzed between a calibration standard with %R between 75-89% (65-79% for Hg; 70-84% for CN) or 111-125% (121-135% for Hg; 116-130% for CN) recovery and nearest good calibration standard. Qualify results <idl %r="" (cn,="" (red-line)="" (uj)="" 65-135%).="" 65-79%).="" 70-130%;="" 70-84%;="" 75-125%="" 75-89%="" as="" ccv="" data="" either="" estimated="" five="" hg,="" icv="" if="" is="" of="" of<="" on="" or="" outside="" qualify="" range="" recovery="" reject="" samples="" side="" td="" the="" unacceptable=""><td>YES NO</td></idl>	YES NO

concentration mode immediately after calibration within $\pm 10\%$ of the true values?

[]

verification standard out of control limits.

A.1.11.	Was continuing calibration performed every 10 samples or every 2 hours?	[]	
	Was ICV for cyanides distilled? []		
AC	TION: If no for any of the above, write in the Contract-Problem/Non-Compliance section of the "Data Assessment Narrative".		
A.1.12	Form II B (CRDL Standards for AA and ICP) -		
A.1.12.	Was a CRDL standard (CRA) analyzed after initial calibration for all AA metals (except Hg)?	[]	
	Was a mid-range calib. verification standard distilled and analyzed for cyanide analysis?	[]	
	Was a 2xCRDL (or 2xIDL when IDL>CRDL) analyzed (CRI) for each ICP run? (Note: CRI for AL, Ba, Ca, Fe, Mg, Na, or K is not required.)	[]	
	ACTION: If no for any of the above, flag as estimated all data falling within the affected ranges. The affected ranges are: AA Analysis - **True Value ± CRDL ICP Analysis - **True Value ± 2CRDL CN Analysis - **True Value ± 0.5 x True Value.		
	value of CRA, CRI or mid-range standard. Substitute IDL for the concentration of the missing mid-range standard from th		
	STANDARD OPERATING PROCEDURE Page	11 of	34
Title:	Evaluation of Metals Data for the Contract Laboratory Program Appendix A.1: Data Assessment - Contract Compliance (Total Review)	Date: Number: Revisio	
A.1.12.2	Was CRI analyzed after ICV/ICB and before the final	YES	NO

	CCV/CCB,	and twice every eigh	t hours of	ICP run?	[]	
	ACTION:	If no, write in Cont Section of the "Data			ice	
A.1.12.3		on each Form IIB all to side the acceptance wi		recoveries tha	it	
	Are CRA	and CRI standards wit	hin contro	l limits:		
			Metals	80 - 120%R?	[]	
	Is mid-r	ange standard within	control lir	mits:		
			Cyanide	80 - 120%R?	[]	
	ACTION:	Flag as estimated al. the affected range is standard is between standard is between standard is between 121-150%; affected range if the reject only positive if the recovery is gothe samples on either the control limits. Flag or reject the fixed data are within the standards are outside	f the recovery flat and reject all a recovery data withing reater than a side of C inal result he affected	very of the ag only positive if the recover leads within the is less than 5 in the affected 150%. Qualifical standard outs only when said ranges and the	y he 0%; range y 50% of tside	
A.1.13	Form III	(Initial and Continua	ing Calibra	ation Blanks)		
A.1.13.1	Present	and complete?			[]	
	For both same an	AA and ICP when both alyte?	are used f	for the	[]	
	Was an i	nitial calibration bla	ank analyze	ed?	[]	
	Was a co every 10 frequent	ntinuing calibration k samples or every 2 ho)?	olank analy ours (which	zed after n ever is more	[]	
		STANDARD OPERATIN	NG PROCEDUR	Æ	Page 12	of 34
Title: Ev	aluation	of Metals Data for the	<u> </u>		Date: Ja	n. 1992

YES NO ACTION: If no, prepare Telephone Record Log, contact laboratory and write in the Contract-Problems/ Non-Compliance section of the "Data Assessment Narrative". A.1.13.2 Circle on each Form III all calibration blank values that are above CRDL (or 2 x IDL when IDL > CRDL). Are all calibration blanks (when IDL<CRDL) less than or equal to the Contract Required Detection Limits (CRDLs)? [] Are all calibration blanks less than two times Instrument Detection Limit (when IDL>CRDL)? [] ACTION: If no for any of the above, flag as estimated (J) positive sample results when raw sample value is less than or equal to calibration blank value analyzed between calibration blank with value over CRDL (or 2xIDL) and nearest good calibration blank. Flag five samples on either side of the calibration blank outside the control limits. A.1.14 FORM III (Preparation Blank) -(Note: The preparation blank for mercury is the same as the calibration blank.) A.1.14.1 Was one prep. blank analyzed for: each Sample Delivery Group (SDG)? each batch of digested samples? each matrix type? both AA and ICP when both are used for the same analyte? ACTION: If no for any of the above, flag as

estimated (J) all the associated positive data <10 x IDLs for which prep. blank

Number:

Revision: 11

HW-2

Contract Laboratory Program

Compliance (Total Review)

Appendix A.1: Data Assessment - Contract

was not analyzed.

NOTE: If only one blank was analyzed for more

than 20 samples, then first 20 samples analyzed

do not have to be flagged as estimated (J).

STANDARD OPERATING PROCEDURE

Title: Evaluation of Metals Data for the

Contract Laboratory Program

Appendix A.1: Data Assessment - Contract

Compliance (Total Review)

Page 13 of 34 Date: Jan. 1992 Number: HW-2 Revision: 11

YES NO A.1.14.2 Is concentration of prep. blank value greater than the CRDL when IDL is less than or equal to CRDL? [] If yes, is the concentration of the sample with the least concentrated analyte less than 10 times the prep.blank? [] ACTION: If yes, reject (red-line) all associated data greater than CRDL concentration but less than ten times the prep. blank value. A.1.14.3 Is concentration of prep. blank value (Form III) less than two times IDL, when IDL is greater than CRDL? ACTION: If no, reject (red-line) all positive sample results when sample raw data are less than 10 times the prep. blank value. A.1.14.4 Is concentration of prep. blank below the negative CRDL? [] ACTION: If yes, reject (red-line) all associated sample results less than 10xCRDL. A.1.15 Form IV (ICP Interference Check Sample) A.1.15.1 Present and complete? [] (NOTE: Not required for furnace AA, flame AA, mercury, cyanide and Ca, Mg, K and Na.) Was ICS analyzed at beginning and end of run (or at least twice every 8 hours)? []

	ACTION: If no, flag as estimated (J) all the samples which AL, Ca, Fe, or Mg is higher than in ICS		
A.1.15.2	Circle all values on each Form IV that are more than \pm 20% of true or established mean value.		
	Are all Interference Check Sample results inside the control limits (\pm 20%)?	[]	
	If no, is concentration of Al, Ca, Fe, or Mg lower than the respective concentration in ICS?	[]	
	STANDARD OPERATING PROCEDURE	Page 14	of 34
	Evaluation of Metals Data for the Contract Laboratory Program Appendix A.1: Data Assessment - Contract Compliance (Total Review)	Date: Ja Number: Revision:	HW-2
	ACTION: If no, flag as estimated (J) those positive results for which ICS recovery is between 121-flag all sample results as estimated if ICS recovery falls within 50-79%; reject (red-line those sample results for which ICS recovery is than 50%; if ICS recovery is above 150%, reject positive results only (not flagged with a "U")	e) s less ct	NQ
A.1.16	Form V A (Spiked Sample Recovery - Pre-Digestion/Pre-Dist (Note: Not required for Ca, Mg, K, and Na (both matrices (soil only.)		
A.1.16.1	Present and complete for: each SDG?	[]	
	each matrix type?	[]	
	each conc. range (i.e. low, med., high)?	[]	
	For both AA and ICP when both are used for the same analyte?	[]	**************************************
	ACTION: If no for any of the above, flag as estimated (J) all the positive data less than four times the spiking levels specified in SOW for which spiked sample was not analyzed.		

 θ_{S, γ_0}

	NOTE:	If one spiked sample was analyzed for mothan 20 samples, then first 20 samples analyzed do not have to be flagged as estimated (J).	ore	
A.1.16.2	Was fie	d blank used for spiked sample?	-	[]
	ACTION:	If yes, flag all positive data less that 4 x spike added as estimated (J) for what field blank was used as spiked sample.	nn nich	
A.1.16.3		on each Form VA all spike recoveries that side control limits (75% to 125%).	:	
	Are all	recoveries within control limits?	[]	
		s sample concentration greater than or e times spike concentration? STANDARD OPERATING PROCEDURE	[]	15 of 34
Co Ap	ntract Lab	f Metals Data for the oratory Program : Data Assessment - Contract Total Review)	Date: Number Revisio	2
			YES	NO
	ACTION:	If yes, disregard spike recoveries for whose concentrations are greater than o to four times spike added. If no, circ analytes on Form V for which sample con is less than four times the spike conce	r equal le those centration	
		ts outside the control limits (75-125%) ith "N" on Form I's and Form VA?	[]	
		If no, write in the Contract - Problem/N Compliance section of "Data Assessment N		
A.1.16.4	Aqueous Are any s	pike recoveries:		
		(a) less than 30%?		[]
		(b) between 30-74%?		[]
		(c) between 126-150%?		[]

		(d)	greater than 150%?			[]
	ACTION:	data; if between 30- aqueous data as esti 126-150%, flag as es aqueous data not fla greater than 150%, r	ject all associated aq 74%, flag all associat mated (J); if between timated (J) all associ gged with a "U"; if eject (red-line) all ata not flagged with a	ed ated		
A.1.16.	5 Soil/Se	diment				
	Are any	spike recoveries:	1.00=			
		(a) less t	han 10%?			[]
		(b) between	n 10-74%?			[]
		(c) between	n 126-200%?			[] _
		(d) greate:	r than 200%?			[]
Title:	Evaluation Contract La Appendix A	ARD OPERATING PROCEDUR of Metals Data for the aboratory Program .1: Data Assessment - (Total Review)	e e		Number:	f 34 Jan. 1992 : HW-2 on: 11
	ACTION:	between 10-74%, flag a if between 126-200%, t data was not flagged w	ect all associated data all associated data as flag as estimated all a with a "U"; if greater data not flagged with	estim associ than	ated 200%,	NO
A.1.17	Form VI	(Lab Duplicates)				
A.1.17.1	l Present	and complete for:	each SDG?		[]	
			each matrix type?		[]	

	each concentration range (i.e. low, med., high)?	[]	
	both AA and ICP when both are used for the same analyte?	[]	
	 ACTION: If no for any the above, flag as estimated (J) all the data ≥CRDL* for which duplicate sample was not analyzed. Note: 1. If one duplicate sample was analyzed for more than 20 samples, then first 20 samples do have to be flagged as estimated. 2. If percent solids for soil sample and its dupli differ by more than 1%, prepare a Form VI for e duplicate pair, report concentrations in ug/L on wet weight basis and calculate RPD or Differ for each analyte. 	cate each	
A.1.17.2	Was field blank used for duplicate analysis?		[]
	ACTION: If yes, flag all data ≥CRDL* as estimated (J) for which field blank was used as duplication	ate.	
A.1.17.3	Are all values within control limits (RPD 20% or difference $\leq \pm \text{CRDL}$)?	[]	
	If no, are all results outside the control limits flagged with an * on Form I's and VI?	[]	
	ACTION: If no, write in the Contract - Problems/Non- Compliance section of "Data Assessment Narras	tive".	
* Subst	itute IDL for CRDL when IDL > CRDL.		
	STANDARD OPERATING PROCEDURE	Page 17	of 34
Con App	luation of Metals Data for the tract Laboratory Program endix A.l: Data Assessment - Contract pliance (Total Review)	Date: Ja Number: Revision	HW-2
	NOTE: 1. RPD is not calculable for an analyte of the sample - duplicate pair when both values are less than IDL. 2. If the result of lab duplicate analyzed by GFAA is rejectable due to coefficient of	YES	NO

correlation of MSA, analytical spike recovery, or duplicate injections criteria, do not apply precision criteria to metals analyzed by GFAA.

A.1.17.4 Aqueous

Circle on each Form VI all values that are:

RPD > 50%, or Difference > CRDL*

Is any RPD greater than 50% where sample and duplicate are both greater than or equal to 5 times *CRDL?

Is any difference** between sample and duplicate greater than *CRDL where sample and/or duplicate is less than 5 times *CRDL?

ACTION: If yes, flag the associated data as estimated.

A.1.17.5 Soil/Sediment

Circle on each Form VI all values that are:

RPD > 100%, or

Difference > 2 x CRDL*

Is any RPD (where sample and duplicate are both greater than or equal to 5 times *CRDL) :

> 100%?

Is any **difference between sample and duplicate (where sample and/or duplicate is less than 5x*CRDL) :

> 2x*CRDL?

STANDARD OPERATING PROCEDURE

Page 18 of 34

[]

Title: Evaluation of Metals Data for the Contract Laboratory Program

Appendix A.1: Data Assessment - Contract

Compliance (Total Review)

Date: Jan. 1992 Number: HW-2

Revision: 11

^{*} Substitute IDL for CRDL when IDL > CRDL.

^{**} Use absolute values of sample and duplicate to calculate the difference.

1

[]

ACTION: If yes, flag the associated data as estimated.

A.1.18 Field Duplicates

A.1.18.1 Were field duplicates analyzed?

ACTION: If yes, prepare a Form VI for each aqueous field duplicate pair. Prepare a Form VI for each soil duplicate pair, if percent solids for sample and its duplicate differ by more than 1%; report concentrations of soils in ug/l on wet weight basis and calculate RPDs or Difference for each analyte.

NOTE: 1. Do not calculate RPD when both values are less than IDL.

Flag all associated data only for field duplicate pair.

A.1.18.2 Aqueous

Circle all values on self prepared Form VI for field duplicates that are:

RPD > 50%, or Difference > CRDL*

Is any RPD greater than 50% where sample and duplicate are both greater than or equal to 5 times *CRDL?

Is any **difference between sample and duplicate greater than *CRDL where sample and/or duplicate is less than 5 times *CRDL?

ACTION: If yes, flag the associated data as estimated.

^{*} Substitute IDL for CRDL when IDL > CRDL.

** [se absolute values of sample and duplicate to calculate the	diffe	renc	e.	
	STANDARD OPERATING PROCEDURE	Page	19	of	3
	Evaluation of Metals Data for the Contract Laboratory Program Appendix A.l: Data Assessment - Contract Compliance (Total Review)	Date: Numbe Revis	r:	HW-	-2
		YES		NO	
A.1.18.3	Soil/Sediment				
	Circle all values on self prepared Form VI for field duplicates that are:				
	RPD >100%, or				
	Difference > 2 x CRDL*				
	Is any RPD (where sample and duplicate are both greater than 5 times *CRDL) : >100%?]]	
	Is any **difference between sample and duplicate (where sample and/or duplicate is less than 5x *CRDL):				
	>2x *CRDL?]]	
	ACTION: If yes, flag the associated data as estimated.				
A.1.19	Form VII (Laboratory Control Sample) (Note: LCS - not required for aqueous Hg and cyanide analyses.)				
A.1.19.1	Was one LCS prepared and analyzed for:				
	each SDG?	[]			
	each batch samples digested/distilled?	[]			
	both AA and ICP when both are used for the same analyte?	[]			
	ACTION: If no for any of the above, prepare Telephone Record Log and contact laboratory for submitta	1			

of results of LCS. Flag as estimated (J) all the data for which LCS was not analyzed.

NOTE: If only one LCS was analyzed for more than 20 samples, then first 20 samples close to LCS do not have to be flagged as estimated.

- * Substitute IDL for CRDL when IDL > CRDL.
- ** Use absolute values of sample and duplicate to calculate the difference.

STANDARD OPERATING PROCEDURE Page 20 of 34 Title: Evaluation of Metals Data for the Date: Jan. 1992 Contract Laboratory Program Number: HW-2 Appendix A.1: Data Assessment - Contract Revision: 11 Compliance (Total Review) YES NO A.1.19.2 Aqueous LCS Circle on each Form VII the LCS percent recoveries outside control limits (80 - 120%) except for aqueous Ag and Sb. Is any LCS recovery: less than 50%? between 50% and 79%? between 121% and 150%? greater than 150%? [] ACTION: Less than 50%, reject (red-line) all data; between 50% and 79%, flag all associated data as estimated (J); between 121% and 150%, flag all positive (not flagged with a "U") results as estimated; greater than 150%, reject all positive results.

A.1.19.3 Solid LCS

NOTE: 1. If "Found" value of LCS is rejectable due to duplicate injections or <u>analytical</u> spike recovery criteria, regardless of LCS recovery, flag the associated data

	2.	as estimated (J). If IDL of an analyte is equal to or greater true value of LCS, disregard the "Action" though LCS is out of control limits.	than below even	
		Is LCS "Found" value higher than the controllimits on Form VII?	1	[]
	ACTION:	If yes, qualify all associated positive dat as estimated.	a	
		Is LCS "Found" value lower than the Control limits on Form VII?		[]
	ACTION:	If yes, qualify all associated data as estimated.		
		STANDARD OPERATING PROCEDURE	Page 21	L of 34
Con App	tract Labora	Metals Data for the Story Program Data Assessment - Contract al Review)	Date: 0 Number: Revision	
			YES	NO
A.1.20	Form IX (IC	P Serial Dilution) -		
	for	al dilution analysis is required only initial concentrations equal to or ter than 10 x IDL.		
A.1.20.1	Was Serial	Dilution analysis performed for: each SDG?	[]	
		each matrix type?	[]	
	each	concentration range (i.e. low, med.)?	[]	
	al 10	no for any of the above, flag as estimated l the positive data $\geq 10 \times IDLs$ or $\geq CRDL$ when $\times IDL \leq CRDL$ for which Serial Dilution Analysis not performed.	n	

Was field blank(s) used for Serial Dilution Analysis? A.1.20.2 **ACTION:** If yes, flag all associated data ≥ 10 x IDL as estimated (J). If $10 \times IDL \leq CRDL$, flag all data ≥ CRDL. A.1.20.3 Are results outside control limit flagged with an "E" on Form I's and Form IX when initial concentration on Form IX is equal to 50 times IDL or greater. [] ACTION: If no, write in the Contract-Problem/Non-Compliance section of the "Data Assessment Narrative". A.1.20.4 Circle on each Form IX all percent difference that are outside the control limits for initial concentrations equal to or greater than 10 x IDLs only. Are any % difference values: > 10%? [] ≥ 100%? STANDARD OPERATING PROCEDURE Page 22 of 34 Title: Evaluation of Metals Data for the Date: Jan. 1992 Contract Laboratory Program Number: HW-2Appendix A.1: Data Assessment - Contract Revision: 11 Compliance (Total Review) YES NO ACTION: Flag as estimated (J) all the associated sample data \geq 10xIDLs (or \geq CRDL when 10xIDL \leq CRDL) for which percent difference is greater than 10%

but less than 100%. Reject (red-line) all the associated sample results equal to or greater than 10xIDLs (or \geq CRDL when $10xIDL \leq$ CRDL) for which PD is greater than or equal to 100%.

Note: Flag or reject on Form I's only the sample results whose associated raw data are \geq 10xIDL (or \geq CRDL when 10xIDL≤ CRDL)

A.1.21	Furnace Atomic Absorbtion (AA) OC Analysis		
A.1.21	.1 Are duplicate injections present in furnace raw data (except during full Method of Standard Addition) for each sample analyzed by GFAA?	[]	
	ACTION: If no, <u>reject</u> the data on Form I's for which duplicate injections were not performed.		
A.1.21.	Do the duplicate injection readings agree within 20% Relative Standard Deviation (RSD) or Coefficient of Variation (CV) for concentration greater than CRDL?	[]	
	Was a dilution analyzed for sample with analytical spike recovery less than 40%?	[]	
	ACTION: If no for any of the above, flag all the associated data as estimated.		
A.1.21.	3 Is *analytical spike recovery outside the control limits (85-115%) for any sample?		[]
	ACTION: If yes, flag as estimated the affected sample r if the recovery is between 10-84%; if the reco between 115-200%, flag the associated positive results as estimated; reject the associated sam results if the recovery is less than 10%; reject positive sample results if the recovery is great than 200%.	very is sample ple t	
* Analy	tical spike is not required on the pre-digestion spiked samp STANDARD OPERATING PROCEDURE		3 of 34
Title:	Evaluation of Metals Data for the Contract Laboratory Program Appendix A.l: Data Assessment - Contract Compliance (Total Review)	Date: Number:	
		YES	NO

NOTE: Reject or flag the data only when the affected sample(s) was not subsequently analyzed by Method of Standard Addition.

A.1.22	Form VIII (Method of Standard Addition Results)	
A.1.22.1	Present? []	
	If no, is any Form I result coded with "S" or a "+"?	[]
	ACTION: If yes, write request on Telephone Record Log and contact laboratory for submittal of Form VIII.	
A.1.22.2	Is coefficient of correlation for MSA less than 0.990 for any sample?	[]
	ACTION: If yes, reject (red-line) the affected data.	
A.1.22.3	Was *MSA required for any sample but not performed?	[]
	Is coefficient of correlation for MSA less than 0.995?	[]
	Are MSA calculations outside the linear range of the calibration curve generated at the beginning of the analytical run?	[]
	ACTION: If yes for any of the above, flag all the associated data as estimated (J).	
A.1.22.4	Was proper quantitation procedure followed correctly as outlined in the SOW on page E-23?	
	ACTION: If no, note exception under Contract Problem/ Non-Compliance section of the "Data Assessment Narrative", and prepare a separate list.	
* MSA	is not required on LCS and prep. blank.	

STANDARD OPERATING PROCEDURE

Page 24 of 34

Date: Jan. 1992

Number: HW-2

Revision: 11

Title: Evaluation of Metals Data for the Contract Laboratory Program Appendix A.1: Data Assessment - Contract Compliance (Total Review)

A.1.22

		YES	NO	
A.1.23	Dissolved/Total or Inorganic/Total Analytes -			
A.1.23.1	Were any analyses performed for dissolved as well as total analytes on the same $sample(s)$.		[]	
	Were any analyses performed for inorganic as well as tot (organic + inorganic) analytes on the same sample(s)?	al	[]	
	NOTE: 1. If yes, prepare a list comparing differences between all dissolved (or inorganic) and total analytes. Compute the differences as a percent of the total analyte only when dissolved concentration is greater than CRDL as well as total concentration. 2. Apply the following questions only if inorganic (or dissolved) results are (i) above CRDL, and (ii) greater than total constituents 3. At least one preparation blank, ICS, and LCS should be analyzed in each analytical run.	-		
A.1.23.2	Is the concentration of any dissolved (or inorganic) analyte greater than its total concentration by more than 10%?		[]	
A.1.23.3	Is the concentration of any dissolved (or inorganic) analyte greater than its total concentration by more than 50%?		[]	-
	ACTION: If more than 10%, flag both dissolved (or inorganic) and total values as estimated (J); if more than 50%, reject (red-line) the data for both values.			
A.1.24	Form I (Field Blank) -			
	(Note: Designate "Field Blank" as such on Form I.)			
A.1.24.1	Circle all field blank values on Form I that are greater than CRDL, (or 2 x IDL when IDL $>$ CRDL).			
	Is field blank concentration less than CRDL (or 2 x IDL when IDL > CRDL) for all parameters			

	of asso	ciated aqueous an STANDARD O	d soil samples? PERATING PROCEDURE		[] Page	25 of 3
	Contract La Appendix A.	of Metals Data fo boratory Program 1: Data Assessme (Total Review)			Date: Number: Revisio	
	due to	other QC criteria			YES	NO
	ACTION:	all associated p than or equal to value. Reject or results that who	except field blank positive sample dat to five times the fine Form I's the soil on converted to ug/than or equal to for value in ug/L.	a less eld blank sample L on wet		
A.1.25	Form X,	XI, XII (Verifica	ation of Instrument	al Parameters	3).	
A.1.25.1	Is verif	ication report pre	esent for:			
		Instrument Dete	ection Limits (quar	terly)?	[]	-
	ICP :	Interelement Corre	ection Factors (ann	ually)?	[]	
		ICP I	Linear Ranges (quar	terly)?	[]	
	ACTION:	If no, contact Th	90 of the lab.			
A.1.25.2		Instrument Detecti for Cyanide.)	on Limits) - (Note	: IDL is not		
A.1.25.2	.1 Are IDLs	present for:	all the ana.	lytes?	[]	
		ć	all the instruments	used?	[]	
	For both analyte?	AA and ICP when b	ooth are used for th	he same	[]	
	ACTION:	If no for any of Telephone Record	the above, prepare Log and contact			

laboratory.

A.1.25.2.2	Is IDL greater than CRDL for any	analyte?			[]
	If yes, is the concentration on analyzed on the instrument whose greater than $5 \times \text{IDL}$.	Form I of the s IDL exceeds CF	sample RDL,	[]	
	STANDARD OPERATING E	ROCEDURE		Page 26	of 34
Co Ap	aluation of Metals Data for the ntract Laboratory Program pendix A.l: Data Assessment - Co mpliance (Total Review)	ntract		Date: Number: Revision	Jan. 1992 HW-2 n: 11
	Action: If no, flag as estimat than five times IDL of IDL exceeds CRDL.	ed all values l the instrument	ess whose	YES	ОИ
A.1.25.3	Form XI (Linear Ranges)				
A.1.25.3.1	Was any sample result higher that of ICP.	n high linear r	ange		[]
	Was any sample result higher tha calibration standard for non-ICP				[] _
	If yes for any of the above, was sample diluted to obtain the res			[]	
	ACTION: If no, flag the result as estimated(J).	reported on For	m I		
A.1.26	Percent Solids of Sediments				
A.1.26.1	Are percent solids in sediment(s): < 50%?		[]	
		< 10%?		[]	

ACTION: If yes, qualify as estimated all the

results of a sample that has per cent solids between 10%-50% (i.e. moisture content between 50%-90%). Reject all the results of a sample that has per cent solids less than 10% (i.e. moisture content

greater than 90%).

NOTE: Reject or flag(J) only the sample results

that were not previously rejected or flaged

due to other QC criteria.

STANDARD OPERATING PROCEDURE

Page 27 of 34

Title: Evaluation of Metals Data for the

Contract Laboratory Program

Appendix A.2: Data Assessment Narrative

Date: Jan. 1992 Number: HW-2

Revision: 11

Case#	Site	 Matrix:	Soil
SDG#	 Lab		Water
Contractor	Reviewer		Other

A.2.1 Validation Flags-

The following flags have been applied in red by the da validator and must be considered by the data user.

J- This flag indicates the result qualified as **estimated**

Red-Line- A red-line drawn through a sample result indicates unu value. The red-lined data are known to contain signi errors based on documented information and must not be

by the data user.

Fully Usable Data-

The results that do not carry "J" or "red-line" are fu

usable.

Contractual Qualifiers- The legend of contractual qualifiers applied by the la
 on Form I's is found on page B-20 of SOW ILM01.0.

A.2.2	The data assessment is given below and on the attached	sheets.
	STANDARD OPERATING PROCEDURE	Page 28 of 34
itle:	Evaluation of Metals Data for the Contract Laboratory Program Appendix A.2: Data Assessment Narrative	Date: Jan. 1992 Number: HW-2 Revision: 11
1.2.2	(continuation)	

_		
_		
_		
-		
_		
_		
	STANDARD OPERATING PROCEDURE	Page 29 of
Title:	Evaluation of Metals Data for the	Date: Jan. 1992
	Contract Laboratory Program Appendix A.2: Data Assessment Narrative	Number: HW-2 Revision: 11
		100 1151011. 11
A.2.2	(continuation)	

e: Evaluation of Metals Data for the	Date: Jan 1992
STANDARD OPERATING PROCEDURE	Page 30 of 34

Title Contract Laboratory Program
Appendix A.2: Data Assessment Narrative

Number: HW-2 Revision: 11

.3 Contract-Probl	em/Non-Compliance	
MMB/ESAT Rviewer:	Signature	Date:
	Signature	

7	Verified by:	Date:			
	STANDARD OPERATING PROCEDURE	Page 31 of 3			
Title:	Evaluation of Metals Data for the Contract Laboratory Program Appendix A.3: Contract Non-Compliance (SMO Report)	Date: Jan. 199 Number: HW-2 Revision: 11			
_	CONTRACT NON-COMPLIANCE (SMO REPORT)				
	Regional Review of Uncontrolled Hazardous Wa Site Contract Laboratory Data Package	aste			
		CASE NO			
Inorganic	opied (laboratory name) data package received at Region II has been reviewed and ce data summarized. The data reviewed included: e No.:	d the quality assura			
Conc. & M	atrix:				
that asso	No.() requires that specific analytical work be dociated reports be provided by the contractor to the Region riteria used to determine the performance were based on a Data Completeness Duplicate Anal Matrix Spike Results Blank Analysis Calibration Standards Results MSA Results	ons, EMSL-LV, and SM an examination of: lysis Results			
Items of	non-compliance with the above contract are described belo	OW.			
Comments:					

Taras		
	Reviewer's Initial Date	
	STANDARD OPERATING PROCEDURE	Page 32 of 34
Title:	Evaluation of Metals Data for the Contract Laboratory Program Appendix A.4: Mailing List for Data Reviewers	Date: Jan. 1992 Number: HW-2 Revision: 11

STANDARD OPERATING PROCEDURE

Title: Evaluation of Metals Data for the

Contract Laboratory Program

Apendix A.5: CLP Data Assessment

Summary Form (Inorganics)

Page 33 of 34

Date: Jan. 1992

Number: HW-2

Revision: 11

STANDARD OPERATING PROCEDURE

Page 34 of 34

Title: Evaluation of Metals Data for the

Contract Laboratory Program

Appendix A.6: CLP Data Assessment Checklist

Date: Jan. 1992 Number: HW-2 Revision: 11

Inorganic Analysis

	INORGANIC	REGIONAL	DATA	ASSESSMENT	Region
CASE NO			_	SITE	
LABORATORY				NO. OF SAMPLES/ MATRIX	
SDG#				REVIEWER (IF NOT ESD)	
SOW#			_	REVIEWER'S NAME	

DPO:	DPO: ACTIONFYI COMPLETION DATE					
		D	ATA ASSESSM	ENT SUMMAR	<u>Y</u>	
10. 11. 12.	HOLDING TIMES CALIBRATIONS BLANKS ICS LCS DUPLICATE ANA MATRIX SPIKE MSA SERIAL DILUTI SAMPLE VERIFI OTHER QC OVERALL ASSES O = Data has no M = Data qualifi Z = Data unaccep X = Problems, bu	LYSIS ON CATION SMENT problems/or qued due to majo table. t do not affec	alified due or problems et data.	e to minor	Hg	
AREA	OF CONCERN:					
NOTAL	BLE PERFORMANCE:					