


National Environmental Methods Index

NEMI Method Summary

Method Number:	504.1
Method Name:	Descriptive Name: EDB,DCBP, and 123TCP in Water Using GCECD Official Name: 1,2-Dibromoethane (EDB), 1,2-Dibromo-3-Chloro-Propane (DCBP), and 1,2,3-Trichloropropane (123TCP) in Water by Microextraction and Gas Chromatography
Revision:	Revision 1.1, 1995
Media:	WATER
Subcategory:	Organic
Source Contact Info:	U.S. EPA Office of Ground Water and Drinking Water Technical Support Center 26 West Martin Luther King Drive Cincinnati, OH 45268-1320
Citation:	Methods for the Determination of Organic Compounds in Drinking Water - Supplement III (EPA/600/R-95-131)  PDF method documentation (80165 Byte file)
Brief Method Summary:	A 35-mL sample is extracted with 2-mL of hexane. The concentrations of analytes in the extract are measured using a gas chromatography (GC) system equipped with a linearized electron capture detector (ECD).
Scope And Application:	This method determines EDB, DCBP, and 123TCP (analytes) in finished drinking water and groundwater.
Applicable Conc Range:	0.03 - 200 ug/L
Interferences:	(A) Solvent contamination: Use high purity solvents, and test solvents to prevent contamination. (B) Extracted interferences: Interference can occur from extracted non-target compounds, with retention times similar to target compounds. (C) DBCM Contamination: DBCM is a common disinfection byproduct which elutes near the EDB retention time window, and can obscure the EDB signal.
QC Requirements:	Initial demonstration of laboratory capability, followed by field reagent blanks (FRB), laboratory reagent blanks (LRB), laboratory fortified blanks (LFB), laboratory fortified sample matrix (LFM), and quality control samples (QCS). A MDL for each analyte must also be determined.
Sample Handling:	Samples are collected in 40 mL bottles and are dechlorinated using sodium thiosulfate. When sampling, the tap is allowed to run for 10 minutes, and collection is done at a flow rate of about 500 mL/min. Samples must be chilled at 4°C upon collection and until analysis, and stored in an area that is free of organic solvent vapors.

Max Holding Time:	14 days; extracts up to 24 hours.
Relative Cost:	\$201 to \$400

Method *EDB,DCBP, and 123TCP in Water Using GCECD (504.1)* has 3 analytes associated with it.

Click on a column header for definition.				
Analyte	Detection Level	Bias	Precision	Spiking Level
1,2-Dibromoethane (106-93-4)	.01 ug/L	117 % Rec (SL)	8 RSD (SL)	.1 ug/L
DBCP (96-12-8)	.01 ug/L	111 % Rec (SL)	1 RSD (SL)	.1 ug/L
1,2,3-Trichloropropane (96-18-4)	.02 ug/L	111 % Rec (SL)	4 RSD (SL)	.1 ug/L

[Close Window](#)