

Location: QA Officer's Office
SOP Files
Wet Chemistry Laboratory

1.0 SCOPE

- 1.1 This SOP is applicable to drinking, surface, and saline waters. This procedure is based on method by HACH which states not approved for EPA reporting.

2.0 SUMMARY OF THE METHOD

- 2.1 This SOP is an spectrophotometric procedure for evaluating Turbidity in liquid samples.

3.0 INTERFERENCES

- 3.1 Presence of floating debris and coarse sediments that settle out rapidly will give low readings. Finely divided air bubbles will affect results in a positive manner.

4.0 APPARATUS AND MATERIALS

- 4.1 100 mL volumetric flask
- 4.2 Hach DR/3000 Spectrophotometer
- 4.3 2 clean cuvettes
- 4.4 Parafilm
- 4.5 Pipettes
- 4.6 Volumetric Flask

5.0 REAGENTS

- 5.1 Deionized (DI) water
- 5.2 Hydrazine Sulfate
- 5.3 Hexamethylenetetramine

6.0 SAMPLE HANDLING AND PRESERVATION

- 6.1 Turbidity is best analyzed as soon as possible from the time the sample is collected. Maximum holding time is 48 hours refrigerated at $4\pm 2^{\circ}\text{C}$.
- 6.2 Samples must be provided in unpreserved, clear plastic bottles.

7.0 PROCEDURE

7.1 Preparation of Standard Solution

- 7.1.1 Stir together 10.0 mL DI Water, 0.500 g Hydrazine Sulfate and 0.050 g Hexamethylenetetramine in a 100 mL volumetric flask until all particles are dissolved in the water.
- 7.1.2 Seal tightly with parafilm and leave it sit in a protected place overnight.
- 7.1.3 Bring standard solution to 100 mL before use.

7.2 Color Development and Measurement

- 7.2.1 Select the stored program for Turbidity by pressing 5 and 9 with the numeric keys and the Stored Program key then press enter.
- 7.2.2 Turn the wavelength dial to 450 nm.
- 7.2.3 Place 25.0 mL of sample in a 25 mL cuvette. If the sample is highly colored, dilute as needed. The upper limit for this analysis is 660 Formazin Turbidity Units (FTU).
- 7.2.4 Prepare a Blank by putting 25.0 mL of DI water in another cuvette.
- 7.2.5 Place the Blank in the spectrophotometer so as the 25 mL mark faces forward. Close the light shield.
- 7.2.6 Press the ZERO and then CONC keys. The spectrophotometer is zeroed for the analysis.
- 7.2.7 Place the sample into the cell holder and close the light shield. Results will be displayed in Formazin turbidity units (FTU).

8.0 QUALITY CONTROL

- 8.1 See Table 1

8.2 Samples are analyzed in batches of 20 or less per QC set. The QC samples that are analyzed per batch are:

- Standard Check
- LCS
- MS
- DUP
- MSD (optional)
- Method Blank

Table 1. Quality Control Requirements (Sample Set = 20 samples)

QC Analysis	Required/ Frequency	Limits	Corrective Action	Corrective Action after Reanalyzing
Method (preparation) Blank	Yes One each set	<MDL or 1/10 Regulatory limit	Remove contamination and rerun	Notify client. Flag data.
Laboratory Control Sample (LCS)	Yes One every 20 samples	90%-110%	Rerun	Notify client. Flag data.
Matrix Duplicate	Yes One each set	RPD<20%	Rerun entire set	Notify client. Flag data.
Matrix Spike	Yes One each set	80%-120%	Analyze by Method of Standard Additions	Notify client. Flag data.
Matrix Spike Duplicate	Yes One every 20 samples Level III	80%-120%	Analyze by Method of Standard Additions	Notify client. Flag data.
Dilution & Rerun	No except if result indicates suppressive interference	Does interference persist?	Yes. Rerun with Method of Standard Additions	Notify client. Flag data.

9.0 DOCUMENTATION

9.1 Turbidity Bench Sheet

- 9.1.1 Analyst
- 9.1.2 Analysis
- 9.1.3 Date Run
- 9.1.4 Detection Limit
- 9.1.5 Wavelength
- 9.1.6 Merit #
- 9.1.7 Dilution
- 9.1.8 Concentration (FTU)
- 9.1.9 Result
- 9.1.10 Spike Concentration
- 9.1.11 % Recovery
- 9.1.12 Batch ID

10.0 METHOD PERFORMANCE

10.1 Precision and accuracy studies are performed on as needed basis. (Ex. new instrument, etc.)

10.2 Method Detection limit studies are performed annually.

11.0 REFERENCES

- 11.1 EPA Water NPDES, Method 180.1, EPA Test Methods, Revision 1982, Turbidity (Nephelometric).
- 11.2 Hach DR/3000 Spectrophotometer Manual.

12.0 SAFETY

- 12.1 Every Laboratory area has eyewash, emergency shower, and fire extinguisher. The metals lab also has dust masks available for use with dust samples.
- 12.2 The air system through out the laboratory area is on a 100% fresh air exchange system, this system exchanges 100% the air in the laboratory area with air from outside 6 times per hour and 30 times per hour when the emergency purge button is hit.
- 12.3 A reference file of material safety data sheets (MSDSs) is available to all personnel.

13.0 WASTE DISPOSAL AND POLLUTION PREVENTION

- 13.1 All laboratory waste must be managed, stored, and disposed in accordance with all federal and state laws and regulations.
- 13.2 Additional information can be found in the Sample Disposal SOP and Merit's Waste Management Plan and Handbook.

14.0 APPROVAL & ISSUE

- 14.1 This section indicates which personnel have read, accepted and approved the SOP. All analysts involved with the SOP must acknowledge their comprehension of the SOP with a signature and a date.

Analyst _____ Date _____

Andy Ball, QA Officer _____ Date _____

Maya V. Murshak, Technical Director _____ Date _____