

34421

Technical Report
for
CBS RECORDS
5152 COLUMBIA DRIVE
CARROLLTON, GA 30117

Chain of Custody Data Required for ETC Data Management Summary Reports					
ETC Sample No.	Company	Facility	Sample Point	Date	Time
J3284	CBS RECORDS	CBSCARGWM	WB9	850723	1030


Denis C. K. Lin, Ph.D.
Vice President
Research and Operations

TABLE OF CONTENTS

Introduction

Table 1: Results and Quality Assurance Data

Table 2: Method Performance Data

Methodology

QA Protocol

Report Appendices

 Appendix A – Mass Spectral Data for Quantitated Compounds

 Appendix A1 – Gas Chromatographic Spectral Data for Quantitated Compounds

 Appendix A2 – HPLC Spectral Data for Quantitated Compounds

 Appendix A3 – Total Aromatics and Derivatized Appendix VIII Compounds by UV/VIS

 Appendix C1 – GC/MS Subsidiary Data – Blank Chromatograms

 Appendix E – Chain of Custody Forms

Introduction

This report contains the analytical results on your water sample, WB9 85/07/23 10:30. Per your request, we analyzed for as many Appendix VIII parameters as presently available methodology permits. USEPA instructs that SW846, July, 1982 "Test Methods for Evaluating Solid Waste" is to be followed in the analysis. However, SW846 does not provide sufficient information to analyze for all the parameters in Appendix VIII. We have performed the analysis based on the methodology described in SW846 and divided the parameters into the following categories:

1. Direct Aqueous Injection for water soluble compounds by GC/MS.
2. Purge and Trap GC/MS method for volatile compounds.
3. Extractable Acid/Base/Neutral/Pesticides compounds by GC/MS.
4. Pesticides by GC/EC.
5. Pesticides by GC/FPD.
6. Herbicides by GC/EC.
7. Polar and thermally unstable compounds by HPLC/UV.
8. Metallic and Organometallic compounds by ICAP, AA, and Cold Vapor AA. (Table 4)
9. Conventional.

Furthermore, there are three categories of compounds that cannot be analyzed directly:

1. Compounds which are unstable in water. (Table A)
2. Not Otherwise Specified (N.O.S.) classes of compounds. (Tables B and D)
3. "Exotic" compounds that require special methods. (Table C)

We analyzed for selected representatives from the N.O.S. classes of compounds.

The Appendix VIII Parameters

As proposed on October 1, 1984, (Federal Register, Vol. 49, Page 38786) there are 13 parameters that are unstable in water; and 10 parameters that are exotic in nature and no satisfactory analytical methods are available. Furthermore, for those parameters that have metals in their compositions it is sufficient to analyze for the metal components only. There are 51 parameters that belong to this category. Using the same rationale 14 parameters can be tested as cyanides. Since one exotic, two unstable, and seven cyanide compounds are also classified as metals, the remaining Appendix VIII parameters that need to be addressed in an analytical scheme total 297.

The analytical scheme designed at ETC for the Appendix VIII parameters follows the rationale suggested by SW-846. Various parameters are categorized according to their chemical and physical characteristics and they are grouped together under analytical methods by which they can be quantified at optional method detection limits. We agree with EPA's approach to the unstable, exotic and metallic compounds. We focus our attention on the remaining 297 parameters.

There are 11 parameters among the 297 that ETC believes should also be classified as "unstable in water" (Table 1). EPA, while acknowledging that there are other unstable compounds, suggests that there is some likelihood they can be found in water. Since EPA insists these compounds should be analyzed, we have decided to incorporate them into our analytical scheme but not setting data acceptance criteria for them.

There are another 33 parameters (Table 2) that ETC believes should be classified as exotic compounds and require special or unusual methods for their analysis. Ten of these compounds are antineoplastic agents or other drugs. Five are uncommon compounds where pertinent information is unavailable. Three are alkaloid poisons or mycotoxins. Two are water soluble dyes. One is a very volatile compound that can explode easily. The remaining 12 compounds, along with the rest of the compounds in this exotic group, possess unusual chemical and physical characteristics which in turn dictate that their analysis at a meaningful detection limit will require individual tailor-made analytical methodology. In most of these cases, EPA recommends to include them in a generalized analytical method without any data to support such recommendations. ETC is obliged to include these parameters into its most

recent Appendix VIII analytical scheme, but strongly urges US EPA to reexamine its position. Perhaps only under unusual circumstances analysis is required for the selected exotic parameters in Table 2.

Among the 29₇ parameters there are 28 that are listed "N.O.S." (not otherwise specified). These 28 parameters, along with aflatoxins, coal tars, creosote, cresols, phthalic acid esters, tetrachloroethane, and toluenediamine, constitute 35 complex mixtures or classes of compounds (Table B). Since some of these parameters represent thousands of individual compounds it is necessary to select representatives for the analysis of these parameters. In as much as possible, we have chosen "priority pollutants" and commonly available reference standard compounds to represent the 35 parameters (see Table B).

Throughout the development of the analytical scheme one of the major obstacles was the availability of reference standard compounds for the 375 parameters.

Results

There are 375 listings in Appendix VIII. To assist you in reviewing the results we have assigned a number to each parameter in the "Grand" Table. Also in the Grand Table we classified the parameters according to their analytical categories and indicated which result table the data and QA data are tabulated on. Depending on the analyses ordered, the quality assurance data may include results from the blank, spiked blank, spiked sample (i.e matrix spike) and replicate sample as well as results from surrogate compound analyses.

In the results table, when a compound or element is present below published MDL it is reported as BMDL (Below Method Detection Limit). When a compound is not present at any detectable concentrations it is reported as ND (Not Detected). Matrix spike and replicate analyses, where included, were performed on samples randomly chosen within each quality assurance batch and are therefore not necessarily spikes and replicates of this report's sample. Surrogate compound recovery data and instrument calibration data are included in the Method Performance Data Tables.

The procedures used in the analysis of the sample are described in this report's methodology section. All analytical procedures within our laboratory are performed within a strictly enforced Quality Assurance Protocol. A description of this Protocol is included in the report. The report appendices include the chain of custody record for your sample and, where appropriate, the gas chromatograms and mass spectra.

TABLE I
PARAMETERS WHICH SHOULD BE UNSTABLE IN WATER

40 Benzotrichloride
53 2-Butanone peroxide
75 1-Chloro-2,3-epoxypropane
79 Chloromethyl methyl ether
138 1,2,3,4-Diepoxybutane
148 Diisopropylfluorophosphate
161 Dimethyl sulfate
220 Maleic anhydride
295 Phorate
298 Phthalic Anhydride
331 Tetraethylthiopyrophosphate

TABLE 2
"EXOTIC" COMPOUNDS THAT REQUIRE SPECIAL METHODS

10 Aflatoxins
15 Mitomycin C
25 Auramine
26 Azaserine
52 Brucine
62 Chlorambucil
86 Citrus red No. 2
98 Cyclophosphamide
99 Daunomycin
103 Diallate
147 3,4-Dihydroxy-alpha-(methylamino)methyl benzyl alcohol
173 2,4-Dithiobiuret
188 Formaldehyde
189 Formic acid
202 Hexaethyl tetraphosphate
241 N-Methyl-N'-nitrosoguanidine
262 4-Nitroquinoline-1-oxide
277 N-Nitrososarcosine
279 Diphosphoramide, octamethyl
281 Endothal
282 Paraldehyde
304 1,3-Propane sultone
306 Propylthiouracil
311 Saccharin and salts
320 Streptozotocin
322 Strychnine and salts
324 2,3,7,8-TCDD
334 Tetranitromethane
343 Thioacetamide
344 Thiosemicarbazide
367 Tri(1-azoidinyl)phosphine sulfide
369 Trypan blue
370 Uracil mustard

GRAND TABLE

COMPOUND NAME	FRACTION	TABLE
1 Acetonitrile	P&T	QR27
2 Acetophenone	A/B/N	QR28
3 Warfarin	HPLC	QR31
4 2-Acetylaminofluorene	A/B/N	QR28
5 Acetyl chloride	UNSTABLE	TBL A
6 1-Acetyl-2-thiourea	HPLC	QR25
7 Acrolein	P&T	QR27
8 Acrylamide	HPLC	QR25
9 Acrylonitrile	P&T	QR27
10 Aflatoxins	CLASS	TBL B
10A Aflatoxins, Total	A/B/N	QR28
11 Aldrin	PEST/HERB	QR24
12 Allyl alcohol	DIR/INJ	QR26
13 Aluminum phosphide	INORG,UNSTABLE	TBL D,A
13A Aluminum	METALS	QR29
14 4-Aminobiphenyl	A/B/N	QR28
15 Mitomycin C	HPLC	QR31
16 5-(Aminomethyl)-3-isoxazolol	A/B/N	QR28
17 Amitrole	HPLC	QR25
18 Aniline	A/B/N	QR28
19 Antimony and Compounds, N.O.S.	CLASS, INORG	TBL B,D
19A Antimony	METALS	QR29
20 Aramite	A/B/N	QR28
21 Arsenic and Compounds, N.O.S.	CLASS, INORG	TBL B,D
21A Arsenic	METALS	QR29
22 Arsenic acid (Orthoarsenic acid)	INORG	TBL D
23 Arsenic pentoxide (Arsenic (V) oxide)	INORG	TBL D
24 Arsenic trioxide (Arsenic (III) oxide)	INORG	TBL D
25 Auramine	A/B/N	QR28
26 Azaserine	HPLC	QR31
27 Barium and Compounds, N.O.S.	CLASS, INORG	TBL B,D
27A Barium	METALS	QR29
28 Barium cyanide	INORG	TBL D
29 Benz[c]acridine	A/B/N	QR28
30 Benz[a]anthracene	A/B/N	QR28
31 Benzene	P&T	QR27
32 Benzenearsonic acid	INORG	TBL D
33 Dichloromethylbenzene	A/B/N	QR28
34 Benzenethiol	A/B/N	QR28
35 Benzidine	HPLC	QR31
36 Benzo[b]fluoranthene	A/B/N	QR28
37 Benzo[j]fluoranthene	A/B/N	QR28
38 Benzo[a]pyrene	A/B/N	QR28
39 p-Benzoquinone	A/B/N	QR28
40 Benzotrichloride	A/B/N	QR28
41 Benzyl chloride	A/B/N	QR28
42 Beryllium and Compounds, N.O.S.	CLASS, INORG	TBL B,D
42A Beryllium	METALS	QR29
43 bis(2-Chloroethoxy)methane	A/B/N	QR28
44 bis(2-Chloroethyl) ether	A/B/N	QR28
45 Chlornaphazine	A/B/N	QR28
46 bis(2-Chloroisopropyl)ether	A/B/N	QR28
47 bis(Chloromethyl)ether	UNSTABLE	TBL A
48 bis(2-ethylhexyl)phthalate	A/B/N	QR28
49 Bromoacetone	P&T	QR27
50 Methyl bromide	P&T	QR27
51 4-Bromophenyl phenyl ether	A/B/N	QR28
52 Brucine	A/B/N	QR28
53 2-Butanone peroxide	A/B/N	QR28

COMPOUND NAME	FRACTION	TABLE
54 Butyl benzyl phthalate	A/B/N	QR28
55 2-sec-Butyl-4,6-dinitrophenol	A/B/N	QR28
56 Cadmium and Compounds, N.O.S.	CLASS, INORG	TBL B,D
56A Cadmium	METALS	QR29
57 Calcium chromate(Chromic acid calcium salt)	INORG	TBL D
57A Calcium	METALS	QR29
58 Calcium cyanide	INORG	TBL D
59 Carbon disulfide	P&T	QR27
60 Carbon oxyfluoride	UNSTABLE	TBL A
61 Chloral	DIR/INJ	QR26
62 Chlorambucil	HPLC	QR31
63 Chlordane	PEST/HERB	QR24
64 Chlorinated Benzenes, N.O.S.	CLASS	TBL B
64A 1,2,3-trichlorobenzene	A/B/N	QR28
64B 2,4,6-trichlorobenzene	A/B/N	QR28
64C 1,2,3,4-tetrachlorobenzene	A/B/N	QR28
64D 1,2,3,5-tetrachlorobenzene	A/B/N	QR28
65 Chlorinated Ethane, N.O.S.	CLASS	TBL B
65A Chloroethane	P&T	QR27
66 Chlorinated Fluorocarbons, N.O.S.	CLASS	TBL B
66A Freon TF	P&T	QR27
67 Chlorinated Naphthalene, N.O.S.	CLASS	TBL B
67A 1-chloronaphthalene	A/B/N	QR28
68 Chlorinated Phenol, N.O.S.	CLASS	TBL B
68A 2,3,5,6-tetrachlorophenol	A/B/N	QR28
68B 2,3,4,5-tetrachlorophenol	A/B/N	QR28
69 Chloroacetaldehyde	DIR/INJ	QR26
70 Chloroalkyl Ethers, N.O.S.	CLASS	TBL B
71 p-Chloroaniline	A/B/N	QR28
72 Chlorobenzene	P&T	QR27
73 Chlorobenzilate	PEST/HERB	QR24
74 p-Chloro-m-cresol	A/B/N	QR28
75 1-Chloro-2,3-epoxypropane	P&T	QR27
76 2-Chloroethyl vinyl ether	P&T	QR27
77 Chloroform	P&T	QR27
78 Methyl chloride	P&T	QR27
79 Chloromethyl methyl ether	P&T	QR27
80 2-Chloronaphthalene	A/B/N	QR28
81 2-Chlorophenol	A/B/N	QR28
82 1-(o-Chlorophenyl)thiourea	HPLC	QR25
83 3-Chloropropionitrile	DIR/INJ	QR26
84 Chromium and Compounds, N.O.S.	CLASS, INORG	TBL B,D
84A Chromium	METALS	QR29
85 Chrysene	A/B/N	QR28
86 Citrus red No. 2	HPLC	QR31
87 Coal Tars	CLASS	TBL B
87A Acenaphthene	A/B/N	QR28
87B Acenaphthalene	A/B/N	QR28
87C Anthracene	A/B/N	QR28
87D Benzo(ghi)perylene	A/B/N	QR28
87E Benzo(k)fluoranthene	A/B/N	QR28
87F Fluorene	A/B/N	QR28
87G Phenanthrene	A/B/N	QR28
87H Pyrene	A/B/N	QR28
88 Copper cyanide	INORG	TBL D
88A Copper	METALS	QR29
89 Creosote	CLASS	TBL B
89A 2-Nitrophenol	A/B/N	QR28
90 Cresols	CLASS	TBL B
90A o-Cresol	A/B/N	QR28
90B m+p-Cresol	A/B/N	QR28

	COMPOUND NAME	FRACTION	TABLE
91	Crotonaldehyde	P&T	QR27
92	Cyanides (soluble salts and complexes) N.O.S.	CLASS, INORG	TBL B,D
92A	Cyanide, Total	METALS	QR29
93	Cyanogen	DIR/INJ	QR26
94	Cyanogen bromide	INORG	TBL D
95	Cyanogen chloride	INORG	TBL D
96	Cycasin	EXOTIC, UNSTABLE	TBL C, A
97	2-Cyclohexyl-4,6-dinitrophenol	A/B/N	QR28
98	Cyclophosphamide	EXOTIC	TBL C
99	Daunomycin	HPLC	QR31
100	4,4'-DDD	PEST/HERB	QR24
101	4,4'-DDE	PEST/HERB	QR24
102	4,4'-DDT	PEST/HERB	QR24
103	Diallate	A/B/N	QR28
104	Dibenz[a,h]acridine	A/B/N	QR28
105	Dibenz[a,j]acridine	A/B/N	QR28
106	Dibenz[a,h]anthracene	A/B/N	QR28
107	7H-Dibenzo[c,g]carbazole	A/B/N	QR28
108	Dibenzo[a,e]pyrene	A/B/N	QR28
109	Dibenzo[a,h]pyrene	A/B/N	QR28
110	Dibenzo[a,i]pyrene	A/B/N	QR28
111	1,2-dibromo-3-chloropropane	P&T	QR27
112	1,2-Dibromoethane	P&T	QR27
113	Dibromomethane	P&T	QR27
114	Di-n-butyl phthalate	A/B/N	QR28
115	1,2-Dichlorobenzene	A/B/N	QR28
116	1,3-Dichlorobenzene	A/B/N	QR28
117	1,4-Dichlorobenzene	A/B/N	QR28
118	Dichlorobenzene, N.O.S.	CLASS	TBL B
119	3,3'-Dichlorobenzidine	HPLC	QR31
120	1,4-Dichloro-2-butene	P&T	QR27
121	Dichlorodifluoromethane	P&T	QR27
122	1,1-Dichloroethane	P&T	QR27
123	1,2-Dichloroethane	P&T	QR27
124	1,2-Trans-dichloroethylene	P&T	QR27
125	Dichloroethylene, N.O.S.	CLASS	TBL B
126	1,1-Dichloroethylene	P&T	QR27
127	Methylene chloride	P&T	QR27
128	2,4-Dichlorophenol	A/B/N	QR28
129	2,6-Dichlorophenol	A/B/N	QR28
130	2,4-D	PEST/HERB	QR24
131	Dichlorophenylarsine	A/B/N	QR28
132	Dichloropropane, N.O.S.	CLASS	TBL B
132A	1,3-Dichloropropane	P&T	QR27
133	1,2-Dichloropropane	P&T	QR27
134	Dichloropropanol, N.O.S.	CLASS	TBL B
134A	Dichloropropanol	DIR/INJ	QR26
135	2,3-Dichloropropene	P&T	QR27
136A	cis-1,3-Dichloropropene	P&T	QR27
136B	trans-1,3-Dichloropropene	P&T	QR27
137	Dieldrin	PEST/HERB	QR24
138	1,2,3,4-Diepoxybutane	P&T	QR27
139	Diethylarsine	P&T	QR27
140	N,N,-Diethylhydrazine	P&T	QR27
141	Carbophenothion	PEST/HERB	QR24
142	O,O-Diethylphosphoric acid, O-p-nitrophenyl ester	A/B/N	QR28
143	Diethyl phthalate	A/B/N	QR28
144	Thionazin	PEST/HERB	QR24
145	Diethylstilbestrol	HPLC	QR25
146	Dihydrosafrole	A/B/N	QR28
147	3,4-Dihydroxy-alpha-(methylamino)methyl benzyl alcohol	HPLC	QR31

	COMPOUND NAME	FRACTION	TABLE
148	Diisopropylfluorophosphate	A/B/N	QR28
149	Dimethoate	PEST/HERB	QR24
150	3,3'-Dimethoxybenzidine	HPLC	QR31
151	p-Dimethylaminoazobenzene	A/B/N	QR28
152	7,12-Dimethylbenz[a]anthracene	A/B/N	QR28
153	3,3'-Dimethylbenzidine	HPLC	QR31
154	Dimethyl carbamoylchloride	UNSTABLE	TBL A
155	1,1-Dimethylhydrazine	DIR/INJ	QR26
156	1,2-Dimethylhydrazine	DIR/INJ	QR26
157	Thiocfanox	A/B/N	QR28
158	alpha-alpha-Dimethylphenethylamine	A/B/N	QR28
159	2,4-Dimethylphenol	A/B/N	QR28
160	Dimethyl phthalate	A/B/N	QR28
161	Dimethyl sulfate	A/B/N	QR28
162	Dinitrobenzene, N.O.S.	CLASS	TBL B
162A	m-Dinitrobenzene	A/B/N	QR28
163	4,6-Dinitro-o-cresol	A/B/N	QR28
164	2,4-Dinitrophenol	A/B/N	QR28
165	2,4-Dinitrotoluene	A/B/N	QR28
166	2,6-Dinitrotoluene	A/B/N	QR28
167	Di-n-octyl phthalate	A/B/N	QR28
168	1,4-Dioxane	DIR/INJ	QR26
169	Diphenylamine	A/B/N	QR28
170	1,2-Diphenylhydrazine	A/B/N	QR28
171	N-Nitrosodi-n-propylamine	A/B/N	QR28
172	Disulfoton	PEST/HERB	QR24
173	2,4-Dithiobiuret	A/B/N	QR28
174A	Endosulfan I	PEST/HERB	QR24
174B	Endosulfan II	PEST/HERB	QR24
175	Endrin	PEST/HERB	QR24
176	Ethyl carbamate	HPLC	QR25
177	Ethyl cyanide	DIR/INJ	QR26
178	Ethylenebisdithiocarbamic acid	EXOTIC	TBL C
179	Ethyleneimine	HPLC	QR25
180	Ethylene oxide	DIR/INJ	QR26
181	Ethylenethiourea	HPLC	QR25
182	Ethyl methacrylate	P&T	QR27
183	Ethyl methanesulfonate	A/B/N	QR28
184	Fluoranthene	A/B/N	QR28
185	Fluorine	UNSTABLE	TBL A
186	2-Fluoroacetamide	UNSTABLE, EXOTIC	TBL A, C
187	Fluoroacetic acid	DIR/INJ	QR26
188	Formaldehyde	P&T	QR27
189	Formic acid	A/B/N	QR28
190	Glycidylaldehyde	DIR/INJ	QR26
191	Halomethane, N.O.S.	CLASS	TBL B
191A	Chlorodibromomethane	P&T	QR27
191B	Dichlorobromomethane	P&T	QR27
192	Heptachlor	PEST/HERB	QR24
193	Heptachlor epoxide	PEST/HERB	QR24
194	Hexachlorobenzene	A/B/N	QR28
195	Hexachlorobutadiene	A/B/N	QR28
196A	Alpha-BHC	PEST/HERB	QR24
196B	Beta-BHC	PEST/HERB	QR24
196C	Gamma-BHC	PEST/HERB	QR24
196D	Delta-BHC	PEST/HERB	QR24
197	Hexachlorocyclopentadiene	A/B/N	QR28
198	Hexachloroethane	A/B/N	QR28
199	Hexachlorohexahydro-endo, endo-dimethanonaphthalene	A/B/N	QR28
200	Hexachlorophene	A/B/N	QR28

COMPOUND NAME	FRACTION	TABLE
201 Hexachloropropene	A/B/N	QR28
202 Hexaethyltetraphosphate	A/B/N	QR28
203 Hydrazine	DIR/INJ	QR26
204 Hydrogen cyanide	SEE 92	TBL B
205 Hydrofluoric acid	UNSTABLE	TBL A
206 Hydrogen sulfide	P&T	QR27
207 Hydroxydimethylarsine oxide	INORG	TBL D
208 Indeno(1,2,3-cd)pyrene	A/B/N	QR28
209 Iodomethane	P&T	QR27
210 Iron dextran	UNSTABLE, EXOTIC	TBL A,C
210A Iron	METAL	QR29
211 Methyl isocyanate	UNSTABLE	TBL A
212 Isobutyl alcohol	DIR/INJ	QR26
213 Isosafrole	A/B/N	QR28
214 Kepone	PEST/HERB	QR24
215 Lasiocarpine	UNSTABLE, EXOTIC	TBL A,C
216 Lead and Compounds, N.O.S.	CLASS, INORG	TBL B,D
216A Lead	METALS	QR29
217 Lead acetate (Acetic acid, lead salt)	INORG	TBL D
218 Lead phosphate (Phosphoric acid, lead salt)	INORG	TBL D
219 Lead subacetate	INORG	TBL D
220 Maleic anhydride	A/B/N	QR28
221 Maleic hydrazide	HPLC	QR25
222 Malononitrile	HPLC	QR25
223 Melphalan	A/B/N	QR28
224 Mercury fulminate	INORG	TBL D
225 Mercury and Compounds, N.O.S.	CLASS, INORG	TBL B,D
225A Mercury	METALS	QR29
226 Methacrylonitrile	DIR/INJ	QR26
227 Methanethiol	DIR/INJ	QR26
228 Methapyrilene	A/B/N	QR28
229 Methomyl	HPLC	QR25
230 Methoxychlor	PEST/HERB	QR24
231 2-Methylaziridine	HPLC	QR25
232 3-Methylcholanthrene	A/B/N	QR28
233 Methyl Chlorocarbonate	UNSTABLE	TBL A
234 4,4'-Methylenebis(2-chloroaniline)	A/B/N	QR28
235 Methyl ethyl ketone	P&T	QR27
236 Methyl hydrazine	DIR/INJ	QR26
237 2-Methylacetonitrile	DIR/INJ	QR26
238 Methyl methacrylate	P&T	QR27
239 Methyl methanesulfonate	A/B/N	QR28
240 Aldicarb	A/B/N	QR28
241 N-Methyl-N'-nitrosoguanidine	A/B/N	QR28
242 Methyl parathion	PEST/HERB	QR24
243 Methylthiouracil	A/B/N	QR28
244 Mustard gas	UNSTABLE, EXOTIC	TBL A,C
245 Naphthalene	A/B/N	QR28
246 1,4-Naphthoquinone	A/B/N	QR28
247 1-Naphthylamine	A/B/N	QR28
248 2-Naphthylamine	A/B/N	QR28
249 1-Naphthyl-2-thiourea	HPLC	QR31
250 Nickel and Compounds, N.O.S.	CLASS, INORG	TBL B,D
250A Nickel	METALS	QR29
251 Nickel carbonyl (Nickel tetracarbonyl)	INORG	TBL D

COMPOUND NAME	FRACTION	TABLE
252 Nickel cyanide (Nickel (II)cyanide)	INORG	TBL D
253 Nicotinic acid	HPLC	QR25
254 Nitric oxide	UNSTABLE, EXOTIC	TBL A,C
255 p-Nitroaniline	A/B/N	QR28
256 Nitrobenzene	A/B/N	QR28
257 Nitrogen dioxide	UNSTABLE	TBL A
258 Nitrogen mustard and hydrochloride salt	EXOTIC, UNSTABLE	TBL C,A
259 Nitrogen mustard N-Oxide and hydrochloride salt	EXOTIC, UNSTABLE	TBL C,A
260 Nitroglycerin	HPLC	QR25
261 4-Nitrophenol	A/B/N	QR28
262 4-Nitroquinoline-1-oxide	A/B/N	QR28
263 Nitrosamines, N.O.S.	CLASS	TBL B
263A N-Nitrosodiphenylamine	A/B/N	QR28
264 N-Nitrosodi-n-butylamine	A/B/N	QR28
265 N-Nitrosodiethanolamine	A/B/N	QR28
266 N-Nitrosodiethylamine	A/B/N	QR28
267 N-Nitrosodimethylamine	A/B/N	QR28
268 N-Nitroso-N-ethylurea	HPLC	QR25
269 N-Nitrosomethylethylamine	A/B/N	QR28
270 N-Nitroso-N-methylurea	HPLC	QR25
271 N-Nitroso-N-methylurethane	A/B/N	QR28
272 N-Nitrosomethylvinylamine	A/B/N	QR28
273 N-Nitrosomorpholine	A/B/N	QR28
274 N-Nitrosonornicotine	A/B/N	QR28
275 N-Nitrosopiperidine	A/B/N	QR28
276 N-Nitrosopyrrolidine	DIR/INJ	QR26
277 N-Nitrososarcosine	A/B/N	QR28
278 5-Nitro-o-toluidine	A/B/N	QR28
279 Octamethylpyrophosphoramido	A/B/N	QR28
280 Osmium tetroxide (Osmium (VIII) oxide)	INORG	TBL D
280A Osmium	METALS	QR29
281 Endothal	A/B/N	QR28
282 Paraldehyde	P&T	QR27
283 Parathion	PEST/HERB	QR24
284 Pentachlorobenzene	A/B/N	QR28
285 Pentachloroethane	P&T	QR27
286 Pentachloronitrobenzene	A/B/N	QR28
287 Pentachlorophenol	A/B/N	QR28
288 Phenacetin	A/B/N	QR28
289 Phenol	A/B/N	QR28
290A m-phenylenediamine	HPLC	QR31
290B o-phenylenediamine	HPLC	QR31
290C p-phenylenediamine	HPLC	QR31
291 Phenylmercury acetate	INORG	TBL D
292 N-Phenylthiourea	HPLC	QR25
293 Phosgene	UNSTABLE	TBL A
294 Phosphine	EXOTIC, UNSTABLE	TBL C, A
295 Phorate	PEST/HERB	QR24
296 Famphur	PEST/HERB	QR24
297 Phthalic acid esters	CLASS	TBL B
298 Phthalic anhydride	A/B/N	QR28
299 2-Picoline	A/B/N	QR28
300 Polychlorinated Biphenyl, N.O.S.	CLASS	TBL B
300A Aroclor 1242	PCB/GC/EC	QR30
300B Aroclor 1254	PCB/GC/EC	QR30
300C Aroclor 1260	PCB/GC/EC	QR30
300D Aroclor 1248	PCB/GC/EC	QR30
300E Aroclor 1232	PCB/GC/EC	QR30
300F Aroclor 1221	PCB/GC/EC	QR30
300G Aroclor 1016	PCB/GC/EC	QR30
301 Potassium cyanide	INORG	TBL D
301A Potassium	METAL	QR29

COMPOUND NAME	FRACTION	TABLE
302 Potassium silver cyanide	INORG	TBL D
303 Pronamide	A/B/N	QR28
304 1,3-Propane sultone	A/B/N	QR28
305 n-Propylamine	DIR/INJ	QR26
306 Propylthiouracil	A/B/N	QR28
307 2-Propyn-1-ol	DIR/INJ	QR26
308 Pyridine	DIR/INJ	QR26
309 Reserpine	HPLC	QR25
310 Resorcinol	A/B/N	QR28
311 Saccharin and salts	A/B/N	QR28
312 Safrole	A/B/N	QR28
313 Selenious acid (Selenium dioxide)	INORG	TBL D
314 Selenium and Compounds, N.O.S.	CLASS, INORG	TBL B,D
314A Selenium	METALS	QR29
315 Selenium sulfide (Sulfur selenide)	INORG	TBL D
316 Selenourea	INORG	TBL D
317 Silver and Compounds, N.O.S.	CLASS, INORG	TBL B,D
317A Silver	METALS	QR29
318 Silver cyanide	INORG	TBL D
319 Sodium cyanide	INORG	TBL D
319A Sodium	METALS	QR29
320 Streptozotocin	HPLC	QR31
321 Strontium sulfide	INORG	TBL D
321A Strontium	METALS	QR29
322 Strychnine and salts	EXOTIC	TBL C
323 1,2,4,5-Tetrachlorobenzene	A/B/N	QR28
324 2,3,7,8-TCDD	EXOTIC	TBL C
325 Tetrachloroethane	CLASS	TBL B
326 1,1,1,2-Tetrachloroethane	P&T	QR27
327 1,1,2,2-Tetrachloroethane	P&T	QR27
328 Tetrachloroethylene	P&T	QR27
329 Carbon tetrachloride	P&T	QR27
330 2,3,4,6-Tetrachlorophenol	A/B/N	QR28
331 Tetraethylthiopyrophosphate	A/B/N	QR28
332 Tetraethyl lead	INORG	TBL D
333 Tetraethylpyrophosphate	PEST/HERB	QR24
334 Tetrinitromethane	P&T	QR27
335 Thallium and Compounds, N.O.S.	CLASS, INORG	TBL B,D
335A Thallium	METALS	QR29
336 Thallic oxide (Thallium (III) oxide)	INORG	TBL D
337 Thallium (I) acetate (Acetic acid, thallium (I) salt)	INORG	TBL D
338 Thallium (I) carbonate	INORG	TBL D
339 Thallium (I) chloride	INORG	TBL D
340 Thallium (I) nitrate	INORG	TBL D
341 Thallium selenite	INORG	TBL D
342 Thallium (I) sulfate	INORG	TBL D
343 Thioacetamide	HPLC	QR31
344 Thiosemicarbazide	HPLC	QR31
345 Thiourea	HPLC	QR25
346 Thiram	A/B/N	QR28
347 Toluene	P&T	QR27
348 Toluenediamine	CLASS	TBL B
348A Toluene-2,4-Diamine	HPLC	QR31
349 O-Toliduidine hydrochloride	A/B/N	QR28
350 Toluene diisocyanate	UNSTABLE	TBL A
351 Toxaphene	PEST/HERB	QR24

COMPOUND NAME	FRACTION	TABLE
352 Bromoform	P&T	QR27
353 1,2,4-Trichlorobenzene	A/B/N	QR28
354 1,1,1-Trichloroethane	P&T	QR27
355 1,1,2-Trichloroethane	P&T	QR27
356 Trichloroethylene	P&T	QR27
357 Trichloromethanethiol	P&T	QR27
358 Trichlorofluoromethane	P&T	QR27
359 2,4,5-Trichlorophenol	A/B/N	QR28
360 2,4,6-Trichlorophenol	A/B/N	QR28
361 2,4,5-T	PEST/HERB	QR24
362 2,4,5-TP (Silvex)	PEST/HERB	QR24
363 Trichloropropane, N.O.S.	CLASS	TBL B
363A 1,1,2-trichloropropane	P&T	QR27
363B 1,2,2-trichloropropane	P&T	QR27
364 1,2,3-Trichloropropane	P&T	QR27
365 0,0,0-Triethyl phosphorothioate	A/B/N	QR28
366 sym-Trinitrobenzene	A/B/N	QR28
367 Tris (1-Azridinyl) phosphine sulfide	A/B/N	QR28
368 Tris(2,3-dibromopropyl) phosphate	A/B/N	QR28
369 Trypan blue	HPLC	QR31
370 Uracil mustard	A/B/N	QR28
371 Vanadic acid, ammonium salt	INORG	TBL D
372 Vanadium pentoxide (Vanadium (V) oxide)	INORG	TBL D
372A Vanadium	METALS	QR29
373 Vinyl chloride	P&T	QR27
374 Zinc cyanide	INORG	TBL D
374A Zinc	METALS	QR29
375 Zinc phosphide	UNSTABLE	TBL A

TABLE A
PARAMETERS WHICH ARE UNSTABLE IN WATER

5 Acetyl chloride
13 Aluminum phosphide
47 bis (Chloromethyl) ether
60 Carbon oxyfluoride
96 Cycasin
154 Dimethyl carbamoylchloride
185 Fluorine
186 2-Fluoroacetamide
205 Hydrofluoric acid
210 Iron dextran
211 Methyl isocyanate
215 Lasiocarpine
233 Methyl chlorocarbonate
244 Mustard gas
254 Nitric oxide
257 Nitrogen dioxide
258 Nitrogen mustard and hydrochloride salt
259 Nitrogen mustard N-Oxide and hydrochloride salt
293 Phosgene
294 Phosphine
350 Toluene diisocyanate
375 Zinc phosphide

TABLE B
COMPOUND CLASSES

CLASS	REPRESENTATIVE	TABLE
19 Antimony and Compounds, N.O.S.	19A Antimony	QR29
21 Arsenic and Compounds, N.O.S.	21A Arsenic	QR29
27 Barium and Compounds, N.O.S.	27A Barium	QR29
42 Beryllium and Compounds, N.O.S.	42A Beryllium	QR29
56 Cadmium and Compounds, N.O.S.	56A Cadmium	QR29
84 Chromium and Compounds, N.O.S.	84A Chromium	QR29
216 Lead and Compounds, N.O.S.	216A Lead	QR29
225 Mercury and Compounds, N.O.S.	225A Mercury	QR29
250 Nickel and Compounds, N.O.S.	250A Nickel	QR29
314 Selenium and Compounds, N.O.S.	314A Selenium	QR29
317 Silver and Compounds, N.O.S.	317A Silver	QR29
335 Thallium and Compounds, N.O.S.	335A Thallium	QR29
10 Aflatoxins	10A Aflatoxins, Total	QR28
64 Chlorinated Benzenes, N.O.S.	64A 1,2,3-trichlorobenzene	QR28
	64B 2,4,6-trichlorobenzene	QR28
	64C 1,2,3,4-tetrachlorobenzene	QR28
	64D 1,2,3,5-tetrachlorobenzene	QR28
	115 1,2-Dichlorobenzene	QR28
	116 1,3-Dichlorobenzene	QR28
	117 1,4-Dichlorobenzene	QR28
	353 1,2,4-Trichlorobenzene	QR28
65 Chlorinated Ethane, N.O.S.	65A Chloroethane	QR27
	354 1,1,1-Trichloroethane	QR27
	355 1,1,2-Trichloroethane	QR27
66 Chlorinated Fluorocarbons, N.O.S.	66A Freon TF	QR27
67 Chlorinated Naphthalene, N.O.S.	67A 1-chloronaphthalene	QR28
68 Chlorinated Phenol, N.O.S.	68A 2,3,5,6-tetrachlorophenol	QR28
	68B 2,3,4,5-tetrachlorophenol	QR28
	128 2,4-Dichlorophenol	QR28
	129 2,6-Dichlorophenol	QR28
	359 2,4,5-Trichlorophenol	QR28
	360 2,4,6-Trichlorophenol	QR28
70 Chloroalkyl Ethers, N.O.S.	44 bis(2-chloroethyl) ether	QR28
	46 bis(2-chloroisopropyl)ether	QR28
	76 2-Chloroethyl vinyl ether	QR27
87 Coal Tars	31 Benzene	QR27
	87A Acenaphthene	QR28
	87B Acenaphthalene	QR28
	87C Anthracene	QR28
	87D 3,4-Benzofluoranthene	QR28
	87E Benzo(ghi)perylene	QR28
	87F Benzo(k)fluoranthene	QR28
	87G Fluorene	QR28
	87H Phenanthrene	QR28
	87I Pyrene	QR28
	90 Cresols	SEE BELOW
	245 Naphthalene	QR28
	289 Phenol	QR28
	347 Toluene	QR27

TABLE B (cont'd)
COMPOUND CLASSES

CLASS	REPRESENTATIVE	TABLE
89 Creosote	74 p-Chloro-m-cresol 89A 2-Nitrophenol 90A o-Cresol 90B m+p Cresol 163 4,6-Dinitro-o-cresol 261 4-Nitrophenol 90A o-Cresol 90B m+p Cresol	QR28 QR28 QR28 QR28 QR28 QR28 QR28 QR28
90 Cresols	92A Cyanide, Total 115 1,2-Dichlorobenzene 116 1,3-Dichlorobenzene 117 1,4-Dichlorobenzene	QR29 QR28 QR28
92 Cyanides, N.O.S. 118 Dichlorobenzene, N.O.S.	126 1,1-Dichloroethylene 132A 1,3-Dichloropropane 133 1,2-Dichloropropane	QR27 QR27 QR27
125 Dichloroethylene, N.O.S. 132 Dichloropropane, N.O.S.	134A Dichloropropanol 162A m-Dinitrobenzene 191A Chlorodibromomethane 191B Dichlorobromoethane	QR26 QR28 QR27 QR27
134 Dichloropropanol, N.O.S. 162 Dinitrobenzene, N.O.S. 191 Halomethane, N.O.S.	263A N-Nitrosodiphenylamine 264 N-Nitrosodi-n-butylamine 265 N-Nitrosodiethanolamine 266 N-Nitrosodiethylamine 267 N-Nitrosodimethylamine 268 N-Nitroso-N-ethylurea 269 N-Nitrosomethylethylamine 270 N-Nitroso-N-methylurea 271 N-Nitroso-N-methylurethane 272 N-Nitrosomethylvinylamine 273 N-Nitrosomorpholine 274 N-Nitrosonornicotine 275 N-Nitrosopiperidine 276 N-Nitrosopyrrolidine 277 N-Nitrososarcosine	QR28 QR28 QR28 QR28 QR28 QR25 QR28 QR25 QR28 QR28 QR28 QR28 QR28 QR28 QR28 QR28
263 Nitrosamines, N.O.S.	48 bis(2-ethylhexyl)phthalate 54 Butyl benzyl phthalate 114 Di-n-Butyl phthalate 143 Diethyl phthalate 160 Dimethyl phthalate 167 Di-n-octyl phthalate	QR28 QR28 QR28 QR28 QR28 QR28
297 Phthalic acid esters	300A Aroclor 1242 300B Aroclor 1254 300C Aroclor 1260 300D Aroclor 1248 300E Aroclor 1232 300F Aroclor 1221 300G Aroclor 1016	QR30 QR30 QR30 QR30 QR30 QR30 QR30
300 Polychlorinated Biphenyl, N.O.S.	326 1,1,1,2-Tetrachloroethane 327 1,1,2,2-Tetrachloroethane 348A Toluene-2,4-Diamine 363A 1,1,2-trichloropropane 363B 1,2,2-Trichloropropane 364 1,2,3-Trichloropropane	QR27 QR27 QR25 QR27 QR27 QR27
325 Tetrachloroethane		
348 Toluenediamine		
363 Trichloropropane, N.O.S.		

TABLE C
"EXOTIC" COMPOUNDS THAT REQUIRE SPECIAL METHODS

96 Cycasin
98 Cyclophosphamide *
178 Ethylenebisdi thiocarbamic acid
186 2-Fluoroacetamide
210 Iron dextran
215 Lasiocarpine
244 Mustard gas
254 Nitric oxide
258 Nitrogen mustard and hydrochloride salt
259 Nitrogen mustard N-Oxide and hydrochloride salt
294 Phosphine
322 Strychnine and salts *
324 2,3,7,8-TCDD **

* No method specified by EPA.

** Method 8280, extra charge required.

TABLE D
INORGANICS
REPRESENTATIVE

TABLE

13 Aluminum phosphide	13A Aluminum	QR29
19 Antimony and Compounds, N.O.S.	19A Antimony	QR29
21 Arsenic and compounds, N.O.S.	21A Arsenic	QR29
22 Arsenic acid (Orthoarsenic acid)	21A Arsenic	QR29
23 Arsenic pentoxide (Arsenic (V) oxide)	21A Arsenic	QR29
24 Arsenic trioxide (Arsenic (III) oxide)	21A Arsenic	QR29
27 Barium and compounds, N.O.S	27A Barium	QR29
28 Barium cyanide	27A Barium, 92A Cyanide, Total	QR29
32 Benzenearsonic acid	21A Arsenic	QR29
42 Beryllium and compounds, N.O.S.	42A Beryllium	QR29
56 Cadmium and compounds, N.O.S.	56A Cadmium	QR29
57 Calcium chromate (Chromic acid, calcium salt)	57A Calcium	QR29
58 Calcium cyanide	84A Chromium 57A Calcium, 92A Cyanide, Total	QR29
84 Chromium and compounds, N.O.S	84A Chromium	QR29
88 Copper cyanide	88A Copper, 92A Cyanide, Total	QR29
92 Cyanides (soluble salts and complexes)	92A Cyanide, Total	QR29
94 Cyanogen bromide	92A Cyanide, Total	QR29
95 Cyanogen chloride	92A Cyanide, Total	QR29
207 Hydroxydimethylarsine oxide	21A Arsenic	QR29
216 Lead and compounds, N.O.S.	216A Lead	QR29
217 Lead acetate (Acetic acid, lead salt)	216A Lead	QR29
218 Lead phosphate (Phosphoric acid, lead salt)	216A Lead	QR29
219 Lead subacetate	216A Lead	QR29
224 Mercury fulminate (Fulminic acid, mercury salt)	225A Mercury	QR29
225 Mercury and compounds, N.O.S.	225A Mercury	QR29
250 Nickel and compounds, N.O.S.	250A Nickel	QR29
251 Nickel carbonyl (Nickel tetracarbonyl)	250A Nickel	QR29
252 Nickel cyanide (Nickel (II)cyanide)	250A Nickel	QR29
280 Osmium tetroxide (Osmium (VIII) oxide)	280A Osmium	QR29
291 Phenylmercury acetate	225A Mercury	QR29
301 Potassium cyanide	301A Potassium 92A Cyanide, Total	QR29
302 Potassium silver cyanide	301A Potassium 317A Silver 92A Cyanide, Total	QR29
313 Selenious acid (Selenium dioxide)	314A Selenium	QR29
314 Selenium and compounds, N.O.S.	314A Selenium	QR29
315 Selenium sulfide (Sulfur selenide)	314A Selenium	QR29
316 Selenourea	314A Selenium	QR29
317 Silver and compounds, N.O.S.	317A Silver	QR29
318 Silver cyanide	317A Silver, 92A Cyanide, Total	QR29
319 Sodium cyanide	319A Sodium, 92A Cyanide, Toatl	QR29
321 Strontium sulfide	321A Strontium	QR29
332 Tetraethyl lead	216A Lead	QR29
335 Thallium and compounds, N.O.S.	335A Thallium	QR29
336 Thallic oxide (Thallium (III) oxide)	335A Thallium	QR29
337 Thallium (I) acetate (Acetic acid, thallium (I) salt)	335A Thallium	QR29
338 Thallium (I) carbonate	335A Thallium	QR29
339 Thallium (I) chloride	335A Thallium	QR29
340 Thallium (I) nitrate	335A Thallium	QR29
341 Thallium selenite	335A Thallium	QR29
342 Thallium (I) sulfate	335A Thallium	QR29
371 Vanadic acid, ammonium salt	371A Vanadium	QR29
372 Vanadium pentoxide (Vanadium (V) oxide)	371A Vanadium	QR29
374 Zinc cyanide	374A Zinc	QR29

J3284

COMMENTS

*Pesticides: Alpha and Gamma BHC results in the QC blank were confirmed on a second column.
4,4-DDD and Endosulfan II are co-eluting compounds.*

AUG 18, 1985

TABLE 1: QUANTITATIVE RESULTS and QUALITY ASSURANCE DATA

RCRA Appendix 8 - Purge & Trap Compounds - GC/MS Analysis Data (QR27)

Chain of Custody Data Required for ETC Data Management Summary Reports										
J3284 CBS RECORDS		CBSCARGMM WB9		850723 1030						
ETC Sample No.	Company	Facility	Sample Point	Date	Time	Elapsed Hours				

Compound *	Results		QC Replicate		QC Blank and Spiked Blank			QC Matrix Spike		
	Sample Concen. ug/l	MDL** ug/l	First ug/l	Second ug/l	Blank Data ug/l	Concen. Added ug/l	% Recov	Unspiked Sample ug/l	Concen. Added ug/l	% Recov
1 Acetonitrile	IND		IND	IND	IND	50	0	IND	50	0
7 Acrolein	ND	100	ND	ND	ND	800	89	ND	800	84
9 Acrylonitrile	ND	100	ND	ND	ND	80	79	ND	80	76
31 Benzene	ND	4.4	ND	ND	ND	18	92	ND	18	104
49 Bromoacetone	-	-	-	-	-	-	-	-	-	-
50 Methyl bromide	ND	10.0	ND	ND	ND	18	97	ND	18	123
59 Carbon disulfide	ND	10	ND	ND	ND	68	22	ND	68	24
65A Chloroethane	ND	10.0	ND	ND	ND	18	80	ND	18	86
66A Freon TF	ND	10.0	ND	ND	ND	68	43	ND	68	51
72 Chlorobenzene	ND	6.0	ND	ND	ND	18	83	ND	18	139
75 1-Chloro-2,3-epoxypropane	-	-	-	-	-	-	-	-	-	-
76 2-Chloroethylvinyl ether	ND	10.0	ND	ND	ND	18	91	ND	18	79
77 Chloroform	ND	1.6	ND	ND	ND	18	85	ND	18	96
78 Methyl chloride	ND	10.0	ND	ND	ND	18	83	ND	18	114
79 Chloromethyl methyl ether	-	-	-	-	-	-	-	-	-	-
91 Crotonaldehyde	ND	10	ND	ND	ND	50	86	ND	50	79
111 1,2-Dibromo-3-chloropropane	ND	10	ND	ND	ND	68	21	ND	68	21
112 1,2-Dibromoethane	ND	10	ND	ND	ND	50	93	ND	50	87
113 Dibromomethane	ND	10	ND	ND	ND	50	89	ND	50	88
120 1,4-Dichloro-2-butene	ND	10	ND	ND	ND	50	90	ND	50	80
121 Dichlorodifluoromethane	ND	10.0	ND	ND	ND	18	79	ND	18	103
122 1,1-Dichloroethane	ND	4.7	ND	ND	ND	18	71	ND	18	107
123 1,2-Dichloroethane	ND	2.8	ND	ND	ND	18	83	ND	18	101
124 1,2-Trans-dichloroethylene	ND	10.0	ND	ND	ND	18	85	ND	18	99
126 1,1-Dichloroethylene	ND	2.8	ND	ND	ND	18	78	ND	18	100
127 Methylene chloride	BMDL	2.8	ND	ND	4.30	18	68	4.20	18	95
132A 1,3-Dichloropropane	ND	10	ND	ND	ND	50	88	ND	50	87
133 1,2-Dichloropropane	ND	6.0	ND	ND	ND	18	83	ND	18	99
135 2,3-Dichloropropene	ND	10	ND	ND	ND	50	89	ND	50	89
136A cis-1,3-Dichloropropylene	ND	10	-	-	-	-	-	-	-	-
136B trans-1,3-Dichloropropylene	ND	10	-	-	-	-	-	-	-	-
138 1,2,3,4-Diepoxybutene	-	-	-	-	-	-	-	-	-	-

*See Ground Table for complete compound name and specifications.

**EPA/ETC established method detection limit.

- = Not detected by current method spectra available.

"ND" = Indeterminate. Standard under spike could not be

detected at current method levels.

Result with neither "MDL" nor "Conc.added" specifies no standard

available. Compound was qualitatively searched for.

AUG 18, 1985

TABLE 1: QUANTITATIVE RESULTS and QUALITY ASSURANCE DATA

RCRA Appendix 8 - Purge & Trap Compounds - GC/MS Analysis Data (QR27)

Chain of Custody Data Required for ETC Data Management Summary Reports										
J3284 CBS RECORDS		CBSCARGWM WB9		850723 1030						
ETC Sample No.	Company	Facility	Sample Point	Date	Time	Elapsed Hours				

Compound *	Results		QC Replicate		QC Blank and Spiked Blank			QC Matrix Spike		
	Sample Concen. ug/l	MDL** ug/l	First ug/l	Second ug/l	Blank Data ug/l	Concen. Added ug/l	% Recov	Unspiked Sample ug/l	Concen. Added ug/l	% Recov
139 Diethylarsine	-	-	-	-	-	-	-	-	-	-
140 N,N-Diethylhydrazine	-	-	-	-	-	-	-	-	-	-
182 Ethyl methacrylate	ND	10	ND	ND	ND	50	89	ND	50	86
188 Formaldehyde	-	-	-	-	-	-	-	-	-	-
191A Chlorodibromomethane	ND	3.1	ND	ND	ND	18	79	ND	18	94
191B Dichlorobromomethane	ND	2.2	ND	ND	ND	18	79	ND	18	98
206 Hydrogen sulfide	-	-	-	-	-	-	-	-	-	-
209 Iodomethane	ND	10	ND	ND	ND	50	91	ND	50	94
235 Methyl ethyl ketone	BMDL	10	ND	ND	ND	140	59	8.85	140	39
238 Methyl methacrylate	ND	10	ND	ND	ND	50	92	ND	50	79
282 Paraldehyde	-	-	-	-	-	-	-	-	-	-
285 Pentachloroethane	ND	10	ND	ND	ND	50	137	ND	50	814
326 1,1,1,2-Tetrachloroethane	ND	10	ND	494	ND	50	32	ND	50	33
327 1,1,2,2-Tetrachloroethane	ND	6.9	ND	ND	ND	18	103	ND	18	110
328 Tetrachloroethylene	ND	4.1	ND	759	ND	18	86	ND	18	44
329 Carbon tetrachloride	ND	2.8	ND	ND	ND	18	89	ND	18	91
334 Tetranitromethane	-	-	-	-	-	-	-	-	-	-
347 Toluene	ND	6.0	ND	1.00	BMDL	18	80	ND	18	107
352 Bromoform	ND	4.7	ND	ND	ND	18	75	ND	18	82
354 1,1,1-Trichloroethane	ND	3.8	ND	ND	BMDL	18	76	ND	18	103
355 1,1,2-Trichloroethane	ND	5.0	ND	ND	ND	18	97	ND	18	111
356 Trichloroethylene	ND	1.9	ND	ND	ND	18	113	3.65	18	48
357 Trichloromethanethiol	-	-	-	-	-	-	-	-	-	-
358 Trichlorofluoromethane	ND	10.0	ND	ND	ND	18	75	ND	18	95
363A 1,1,2-Trichloropropane	ND	10	ND	ND	ND	50	32	ND	50	33
363B 1,2,2-Trichloropropane	ND	-	ND	ND	ND	0	-	ND	0	-
364 1,2,3-Trichloropropane	ND	10	ND	ND	ND	140	48	1.43	140	40
373 Vinyl chloride	ND	10.0	ND	ND	ND	18	83	ND	18	109

*See Grund Table for complete compound name and specifications.
 **EPA/ETC established method detection limit.
 "ND" specifies no standard or reference spectra available.
 "BMDL" is indeterminate. Standard and/or spikes could not be determined due to low signal/noise ratio.
 Result with neither "MDL" nor "Conc. Added" specifies no standard available, compound was qualitatively searched for.

SEP 5, 1985

TABLE 1: QUANTITATIVE RESULTS and QUALITY ASSURANCE DATA

RCRA Appendix 8 - Acid/B/N/Pest Compounds - GC/MS Analysis Data (QR28)

Chain of Custody Data Required for ETC Data Management Summary Reports									
J3284 CBS RECORDS		CBSCARGWM WB9		850723 1030					
ETC Sample No.	Company	Facility	Sample Point	Date	Time	Elapsed Hours			

Compound *	Results		QC Replicate		QC Blank and Spiked Blank			QC Matrix Spike		
	Sample Concen. ug/l	MDL** ug/l	First ug/l	Second ug/l	Blank Data ug/l	Concen. Added ug/l	% Recov	Unspiked Sample ug/l	Concen. Added ug/l	% Recov
2 Acetophenone	ND	10	ND	ND	ND	100	95	ND	101	129
4 2-Acetylaminofluorene	ND	10	ND	ND	ND	300	141	ND	303	99
10A Aflatoxins, Total	-	-	-	-	-	-	-	-	-	-
14 4-Aminobiphenyl	ND	10	ND	ND	ND	100	91	ND	101	110
16 5-(Aminomethyl)-3-isoxazolo	IND	-	IND	IND	IND	300	0	ND	303	0
18 Aniline	ND	10	ND	ND	ND	100	92	ND	101	107
20 Aramite	-	-	-	-	-	-	-	-	-	-
25 Auramine	-	-	-	-	-	-	-	-	-	-
29 Benz[c]acridine	-	-	-	-	-	-	-	-	-	-
30 Benz[a]anthracene	ND	10	ND	ND	ND	100	91	ND	101	153
33 Dichloromethylbenzene	ND	10	ND	ND	ND	100	7	ND	101	42
34 Benzenethiol	ND	10	ND	ND	ND	300	25	ND	303	25
36 Benzo[b]fluoranthene	ND	10	ND	ND	ND	100	132	ND	101	95
37 Benzo[j]fluoranthene	-	-	-	-	-	-	-	-	-	-
38 Benzo[a]pyrene	ND	2.5	ND	ND	ND	100	126	ND	101	108
39 p-Benzoquinone	IND	-	IND	IND	IND	100	0	IND	101	24
40 Benzotrichloride	-	-	-	-	-	-	-	-	-	-
41 Benzyl chloride	ND	10	ND	ND	ND	100	77	ND	101	104
43 bis(2-Chloroethoxy)methane	ND	5.4	ND	ND	ND	100	82	ND	101	117
44 bis(2-Chloroethyl) ether	ND	5.8	ND	ND	ND	100	72	ND	101	143
45 Chlornaphazine	-	-	-	-	-	-	-	-	-	-
46 bis(2-Chloroisopropyl)ether	ND	5.8	ND	ND	ND	100	88	ND	101	118
48 bis(2-Ethylhexyl)phthalate	ND	10	1.18	801	ND	100	98	3.68	101	159
51 4-Bromophenyl phenyl ether	ND	1.9	ND	ND	ND	100	115	ND	101	154
52 Brucine	-	-	-	-	-	-	-	-	-	-
53 2-Butanone peroxide	-	-	-	-	-	-	-	-	-	-
54 Butyl benzyl phthalate	ND	10	ND	ND	ND	100	91	ND	101	148
55 2-sec-Butyl-4,6-dinitrophen	ND	10	ND	ND	ND	300	96	ND	301	139
64A 1,2,3-Trichlorobenzene	ND	10	ND	ND	ND	100	88	ND	101	95
64B 1,3,5-Trichlorobenzene	ND	10	ND	ND	ND	100	92	ND	101	97
64C 1,2,3,4-Tetrachlorobenzene	ND	10	ND	ND	ND	100	111	ND	101	116
64D 1,2,3,5-Tetrachlorobenzene	ND	10	ND	ND	ND	100	98	ND	101	105

*See Grand Table for complete compound name and specifications.

**EPR/ETC established method detection limit.

***No standard or reference spectra available.

****ND = Not detected. Standard and spikes could not be

detected at current method levels.

Result with neither "MDL" nor "Conc.added" specifies no standard

available; compound was qualitatively searched for.

SEP 5, 1985

TABLE 1: QUANTITATIVE RESULTS and QUALITY ASSURANCE DATA

RCRA Appendix 8 – Acid/B/N/Pest Compounds – GC/MS Analysis Data (QR28)

Chain of Custody Data Required for ETC Data Management Summary Reports									
J3284 CBS RECORDS		CBSCARGWM W89		850723 1030					
ETC Sample No.	Company	Facility	Sample Point	Date	Time	Elapsed	Hours		

Compound *	Results			QC Replicate		QC Blank and Spiked Blank			QC Matrix Spike		
	Sample Concen. ug/l	MDL** ug/l		First ug/l	Second ug/l	Blank Data ug/l	Concen. Added ug/l	% Recov	Unspiked Sample ug/l	Concen. Added ug/l	% Recov
67A 1-chloronaphthalene	ND	10		ND	ND	ND	100	317	ND	101	374
68A 2,3,5,6-Tetrachlorophenol	ND	10		ND	ND	ND	100	112	ND	101	56
68B 2,3,4,5-Tetrachlorophenol	ND	10		ND	ND	ND	100	122	ND	101	107
71 p-Chloroaniline	ND	10		ND	ND	ND	300	31	ND	301	44
74 p-Chloro-m-cresol	ND	3.0		ND	ND	ND	100	65	ND	101	85
80 2-Chloronaphthalene	ND	1.9		ND	ND	ND	100	54	ND	101	52
81 2-Chlorophenol	ND	3.3		ND	ND	ND	100	74	ND	101	106
85 Chrysene	ND	2.5		ND	ND	ND	100	90	ND	101	121
87A Acenaphthene	ND	1.9		ND	ND	ND	100	92	ND	101	118
87B Acenaphthylene	ND	3.5		ND	ND	ND	100	98	ND	101	123
87C Anthracene	ND	1.9		ND	ND	ND	100	109	ND	101	149
87D Benzo(ghi)perylene	ND	4.1		ND	ND	ND	0	-	ND	0	-
87E Benzo(k)fluoranthene	ND	3.5		ND	ND	ND	100	102	ND	101	116
87F Fluorene	ND	1.9		ND	ND	ND	100	88	ND	101	92
87G Phenanthrene	ND	5.5		ND	ND	ND	100	114	ND	101	158
87H Pyrene	ND	1.9		ND	ND	ND	100	87	ND	101	152
89A 2-Nitrophenol	ND	3.6		ND	ND	ND	100	86	ND	101	104
90A o-Cresol	ND	10		ND	ND	ND	100	101	ND	101	88
90B m+p-Cresol	ND	10		ND	ND	ND	200	61	ND	202	65
97 2-Cyclohexyl-4,6-dinitrophe	-	-	-	-	-	-	-	-	-	-	-
103 Diallate	-	-	-	-	-	-	-	-	-	-	-
104 Dibenz[a,h]acridine	-	-	-	-	-	-	-	-	-	-	-
105 Dibenz[a,j]acridine	ND	10		ND	ND	ND	100	113	ND	101	58
106 Dibenz[a,h]anthracene	ND	10		ND	ND	ND	0	-	ND	0	-
107 7H-Dibenzo[c,g]carbazole	-	-	-	-	-	-	-	-	-	-	-
108 Dibenzo[a,e]pyrene	-	-	-	-	-	-	-	-	-	-	-
109 Dibenzo[a,h]pyrene	-	-	-	-	-	-	-	-	-	-	-
110 Dibenzo[a,i]pyrene	-	-	-	-	-	-	-	-	-	-	-
114 Di-n-butyl phthalate	ND	10	2.30	1.46	BMDL	ND	100	105	1	101	0
115 1,2-Dichlorobenzene	ND	1.9	ND	ND	ND	ND	100	63	ND	101	83
116 1,3-Dichlorobenzene	ND	1.9	ND	ND	ND	ND	100	50	ND	101	72
117 1,4-Dichlorobenzene	ND	4.4	ND	ND	ND	ND	100	62	ND	101	75

*See Grand Table for complete compound name and specifications.
**EPA/ETC established method detection limit.
***specifies no standard or reference spectra available.
****=indeterminate. Standard and/or spikes could not be detected at current method levels.
Result with asterisk indicates no standard was available, compound was qualitatively searched for.

SEP 5, 1985

TABLE 1: QUANTITATIVE RESULTS and QUALITY ASSURANCE DATA

RCRA Appendix 8 - Acid/B/N/Pest Compounds - GC/MS Analysis Data (QR28)

Chain of Custody Data Required for ETC Data Management Summary Reports										
J3284 CBS RECORDS		CBSCARGWM WB9		850723 1030						
ETC Sample No.	Company	Facility	Sample Point	Date	Time	Elapsed Hours				

Compound *	Results			QC Replicate		QC Blank and Spiked Blank			QC Matrix Spike		
	Sample Concen. ug/l	MDL** ug/l	First ug/l	Second ug/l	Blank Data ug/l	Concen. Added ug/l	% Recov	Unspiked Sample ug/l	Concen. Added ug/l	% Recov	
128 2,4-Dichlorophenol	ND	2.7	ND	ND	ND	100	81	ND	101	95	
129 2,6-Dichlorophenol	ND	10	ND	ND	ND	100	89	ND	101	105	
131 Dichlorophenylarsine	-	-	-	-	-	-	-	-	-	-	
142 0,0-Diethylphosphoric acid	-	-	-	-	-	-	-	-	-	-	
143 Diethyl phthalate	ND	10	ND	ND	ND	100	91	ND	101	99	
146 Dihydrosafrole	-	-	-	-	-	-	-	-	-	-	
148 Diisopropylfluorophosphate	-	-	-	-	-	-	-	-	-	-	
151 p-Dimethylaminoazobenzene	ND	10	ND	ND	ND	100	114	ND	101	141	
152 7,12-Dimethylbenz[a]anthrac	ND	10	ND	ND	ND	100	87	ND	101	68	
157 Thiophanox	IND		IND	IND	IND	300	0	IND	303	0	
158 a-a-Dimethylphenethylamine	IND		IND	IND	IND	100	0	IND	101	0	
159 2,4-Dimethylphenol	ND	2.7	ND	ND	ND	100	90	ND	101	89	
160 Dimethyl phthalate	ND	10	ND	ND	ND	100	92	ND	101	105	
161 Dimethyl sulfate	-	-	-	-	-	-	-	-	-	-	
162A m-Dinitrobenzene	ND	10	ND	ND	ND	100	93	ND	101	109	
163 4,6-Dinitro-o-cresol	ND	24	ND	ND	ND	100	54	ND	101	81	
164 2,4-Dinitrophenol	ND	42	ND	ND	ND	100	52	ND	101	89	
165 2,4-Dinitrotoluene	ND	5.8	ND	ND	ND	100	75	ND	101	83	
166 2,6-Dinitrotoluene	ND	1.9	ND	ND	ND	100	95	ND	101	113	
167 Di-n-octyl phthalate	ND	10	ND	ND	ND	100	109	ND	101	146	
169 Diphenylamine	ND	10	ND	ND	ND	100	0	ND	101	0	
170 1,2-Diphenylhydrazine	ND	10	ND	ND	ND	100	82	ND	101	96	
171 N-Nitrosodi-n-propylamine	ND	10	ND	ND	ND	100	29	ND	101	127	
173 2,4-Dithiobiuret	-	-	-	-	-	-	-	-	-	-	
183 Ethyl methanesulfonate	ND	10	ND	ND	ND	100	79	ND	101	73	
184 Fluoranthene	ND	2.2	ND	ND	ND	100	90	ND	101	118	
189 Formic acid	-	-	-	-	-	-	-	-	-	-	
194 Hexachlorobenzene	ND	1.9	ND	ND	ND	100	97	ND	101	132	
195 Hexachlorobutadiene	ND	.91	ND	ND	ND	100	55	ND	101	72	
197 Hexachlorocyclopentadiene	ND	10	ND	ND	ND	0	-	ND	0	-	
198 Hexachloroethane	ND	1.6	ND	ND	ND	100	56	ND	101	58	
199 Hexachlorohexahydro-endo, e	-	-	-	-	-	-	-	-	-	-	

*See Grand Table for complete compound name and specifications.

**EPA/ETC established method detection limit.

***specifies no standard or reference spectra available.

****Indicates indeterminate. Standard and/or spikes could not be

Result with neither "MDL" nor "Conc.added" specifies no standard available, compound was qualitatively searched for.

SEP 5, 1985

TABLE 1: QUANTITATIVE RESULTS and QUALITY ASSURANCE DATA

RCRA Appendix 8 – Acid/B/N/Pest Compounds – GC/MS Analysis Data (QR 28)

Chain of Custody Data Required for ETC Data Management Summary Reports									
J3284 CBS RECORDS		CBSCARGWM WB9		850723 1030					
ETC Sample No.	Company	Facility	Sample Point	Date	Time	Elapsed Hours			

Compound *	Results			QC Replicate		QC Blank and Spiked Blank			QC Matrix Spike		
	Sample Concen. ug/l	MDL** ug/l	First ug/l	Second ug/l	Blank Data ug/l	Concen. Added ug/l	% Recov	Unspiked Sample ug/l	Concen. Added ug/l	% Recov	
200 Hexachlorophene	ND	10	ND	ND	ND	100	117	ND	101	145	
201 Hexachloropropene	ND	10	ND	ND	ND	300	52	ND	303	73	
202 Hexaethyltetraphosphate	-	-	-	-	-	-	-	-	-	-	
208 Indeno(1,2,3-cd)pyrene	ND	10	ND	ND	ND	0	-	ND	0	-	
213 Isosafrole	ND	10	ND	ND	ND	300	184	ND	303	217	
220 Maleic anhydride	-	-	-	-	-	-	-	-	-	-	
223 Melphalan	IND	-	IND	IND	IND	300	0	IND	303	0	
228 Methapyrilene	-	-	-	-	-	-	-	-	-	-	
232 3-Methylcholanthrene	ND	10	ND	ND	ND	100	103	ND	101	56	
234 4,4'-Methylenebis(2-chloro)	ND	10	ND	ND	ND	300	127	ND	303	65	
239 Methyl methanesulfonate	ND	10	ND	ND	ND	100	13	ND	101	97	
240 Aldicarb	ND	10	ND	ND	ND	100	12	ND	101	162	
241 N-Methyl-N'-nitrosoguanidin	-	-	-	-	-	-	-	-	-	-	
243 Methylthiouracil	-	-	-	-	-	-	-	-	-	-	
245 Naphthalene	ND	1.6	ND	ND	ND	100	76	ND	101	101	
246 1,4-Naphthoquinone	ND	10	ND	ND	ND	100	37	ND	101	47	
247 1-Naphthylamine	ND	10	ND	ND	ND	100	50	ND	101	34	
248 2-Naphthylamine	ND	10	ND	ND	ND	100	36	ND	101	17	
255 p-Nitroaniline	ND	10	ND	ND	ND	300	115	ND	303	178	
256 Nitrobenzene	ND	1.9	ND	ND	ND	100	88	ND	101	113	
261 4-Nitrophenol	ND	2.4	ND	ND	ND	100	54	ND	101	41	
262 4-Nitroquinoline-1-oxide	-	-	-	-	-	-	-	-	-	-	
263A N-Nitrosodiphenylamine	ND	1.9	ND	ND	ND	100	176	ND	101	155	
264 N-Nitrosodi-n-butylamine	ND	10	ND	ND	ND	100	87	ND	101	103	
265 N-Nitrosodiethanolamine	IND	-	IND	IND	IND	300	0	IND	303	0	
266 N-Nitrosodiethylamine	ND	10	ND	ND	ND	300	76	ND	303	99	
267 N-Nitrosodimethylamine	ND	10	ND	ND	ND	0	-	ND	0	-	
269 N-Nitrosomethyl ethylamine	IND	10	IND	IND	IND	300	0	IND	303	0	
271 N-Nitroso-N-methylurethane	ND	10	ND	ND	ND	300	43	ND	303	84	
272 N-Nitrosomethyl vinylamine	-	-	-	-	-	-	-	-	-	-	
273 N-Nitrosomorpholine	ND	10	ND	ND	ND	300	55	ND	303	55	
274 N-Nitrosonornicotine	-	-	-	-	-	-	-	-	-	-	

*See Grand Table for complete compound name and specifications.

**EPA/ETC established method detection limit.

***First five significant figures available.

****Indeterminate. Standard and/or spikes could not be detected at current method levels.

Result with neither "MDL" nor "Conc.added" specifies no standard available; compound was qualitatively searched for.

SEP 5, 1985

TABLE 1: QUANTITATIVE RESULTS and QUALITY ASSURANCE DATA

RCRA Appendix 8 - Acid/B/N/Pest Compounds - GC/MS Analysis Data (QR28)

Chain of Custody Data Required for ETC Data Management Summary Reports											
J3284 CBS RECORDS			CBSCARGWM WB9			850723 1030					
ETC Sample No.	Company	Facility	Sample Point	Date	Time	Elapsed Hours					
Compound *	Results	QC Replicate	QC Blank and Spiked Blank	QC Matrix Spike							
	Sample Concen. ug/l	MDL** ug/l	First ug/l	Second ug/l	Blank Data ug/l	Concen. Added ug/l	% Recov	Unspiked Sample ug/l	Concen. Added ug/l	% Recov	
275 N-Nitrosopiperidine	ND	10	ND	ND	ND	100	93	ND	101	143	
277 N-Nitrososarcosine	-	-	-	-	-	-	-	ND	-	-	
278 5-Nitro-o-toluidine	ND	10	ND	ND	ND	300	102	ND	303	163	
279 Octamethylpyrophosphoramide	-	-	-	-	-	-	-	-	-	-	
281 Endothal	-	-	-	-	-	-	-	-	-	-	
284 Pentachlorobenzene	ND	10	ND	ND	ND	100	93	ND	101	116	
286 Pentachloronitrobenzene	ND	10	ND	ND	ND	300	82	ND	303	110	
287 Pentachlorophenol	ND	3.6	ND	ND	ND	100	149	ND	101	164	
288 Phenacetin	ND	10	ND	ND	ND	100	113	ND	101	164	
289 Phenol	ND	1.5	ND	ND	ND	100	84	ND	101	81	
298 Phthalic anhydride	-	-	-	-	-	-	-	-	-	-	
299 2-Picoline	ND	10	ND	ND	ND	100	64	ND	101	101	
303 Pronamide	ND	10	ND	ND	ND	100	99	ND	101	130	
304 1,3-Propane sultone	-	-	-	-	-	-	-	-	-	-	
306 Propylthiouracil	-	-	-	-	-	-	-	-	-	-	
310 Resorcinol	IND	-	IND	IND	IND	300	0	IND	303	0	
311 Saccharin and salts	-	-	-	-	-	-	-	-	-	-	
312 Safrole	ND	10	ND	ND	ND	300	144	ND	303	114	
323 1,2,4,5-Tetrachlorobenzene	ND	10	ND	ND	ND	100	203	ND	101	92	
330 2,3,4,6-Tetrachlorophenol	ND	10	ND	ND	ND	300	65	ND	303	138	
331 Tetraethylidithiopyrophospha	-	-	-	-	-	-	-	-	-	-	
346 Thiram	-	-	-	-	-	-	-	-	-	-	
349 O-Toluidine	IND	-	ND	ND	ND	300	64	ND	303	69	
353 1,2,4-Trichlorobenzene	ND	1.9	ND	ND	ND	100	71	ND	101	90	
359 2,4,5-Trichlorophenol	ND	10	ND	ND	ND	100	91	ND	101	120	
360 2,4,6-Trichlorophenol	ND	2.7	ND	ND	ND	100	71	ND	101	96	
365 0,0,0-Triethyl phosphorothio	-	-	-	-	-	-	-	-	-	-	
366 sym-Trinitrobenzene	-	-	-	-	-	-	-	-	-	-	
367 Tris(1-Azridinyl)phosphine	-	-	-	-	-	-	-	-	-	-	
368 Tris{2,3-dibromopropyl}phos	IND	-	IND	IND	IND	300	0	IND	303	0	
370 Uracil mustard	-	-	-	-	-	-	-	-	-	-	

*See Grand Table for complete compound name and specifications.

**EPA/ETC established method detection limit.

"-":specifies no standard or reference spectra available.

"IND":Indeterminate. Standard and matrix spikes could not be

detected during method level search.

Result with neither "MDL" nor "Conc.added" specifies no standard available; compound was qualitatively searched for.

AUG 27, 1985

TABLE 1: QUANTITATIVE RESULTS and QUALITY ASSURANCE DATA
BASE/NEUTRAL COMPOUNDS – GC/MS ANALYSIS DATA (QR03)

Chain of Custody Data Required for ETC Data Management Summary Reports												
J3284 CBS RECORDS		CBSCARGMM W89		850723 1030				Facility	Sample Point	Date	Time	Elapsed Hours
ETC Sample No.	Company											

NPDES Number	Compound	Results		QC Replicate		QC Blank and Spiked Blank			QC Matrix Spike		
		Sample Concen. ug/l	MDL ug/l	First ug/l	Second ug/l	Blank Data ug/l	Concen. Added ug/l	% Recov.	Unspiked Sample ug/l	Concen. Added ug/l	% Recov.
	Formaldehyde	ND	1000	-	-	ND	-	-	-	-	-

A Standard QA procedures not applicable to these compounds.

AUG 18, 1985

TABLE 1: QUANTITATIVE RESULTS and QUALITY ASSURANCE DATA
RCRA Appendix 8 - Pest & Herb Compounds - GC Analysis Data (QR24)

Chain of Custody Data Required for ETC Data Management Summary Reports									
J3284 CBS RECORDS			CBSCARGMM WB9		850723 1030				
ETC Sample No.	Company	Facility	Sample Point	Date	Time	Elapsed Hours			

Compound *	Results		QC Replicate		QC Blank and Spiked Blank			QC Matrix Spike		
	Sample Concen. ug/l	MDL** ug/l	First ug/l	Second ug/l	Blank Data ug/l	Concen. Added ug/l	% Recov	Unspiked Sample ug/l	Concen. Added ug/l	% Recov
11 Aldrin	ND	.10	ND	ND	ND	0.50	70	ND	0.56	98
63 Chlordane	ND	1.0	ND	ND	ND	5.0	136	ND	5.6	130
73 Chlorobenzilate	-	-	-	-	-	-	-	-	-	-
100 4,4'-DDD	ND	.10	ND	ND	ND	0.50	0	ND	0.56	127
101 4,4'-DDE	ND	.10	ND	ND	ND	0.50	98	ND	0.56	100
102 4,4'-DDT	ND	.10	ND	ND	ND	0.50	0	ND	0.56	108
130 2,4-D	ND	50	ND	ND	ND	100	106	ND	100	103
137 Dieldrin	ND	.10	ND	ND	ND	0.50	118	ND	0.56	107
141 Carbophenothion	-	-	-	-	-	-	-	-	-	-
144 Thionazin	ND	1.0	ND	ND	ND	5.0	84	ND	5.6	93
149 Dimethoate	ND	5.2	ND	ND	ND	5.0	74	ND	5.6	70
172 Disulfoton	ND	5.2	ND	ND	ND	25	86	ND	28	90
174A Endosulfan I	ND	.10	ND	ND	ND	0.50	0	ND	0.56	98
174B Endosulfan II	ND	.10	ND	ND	ND	0.50	0	ND	0.56	127
175 Endrin	ND	.10	ND	ND	ND	0.50	0	ND	0.56	131
192 Heptachlor	ND	.10	ND	ND	ND	0.50	134	ND	0.56	101
193 Heptachlor epoxide	ND	.10	ND	ND	ND	0.50	116	ND	0.56	109
196A Alpha-BHC	ND	.10	ND	ND	ND	0.50	86	ND	0.56	96
196B Beta-BHC	ND	.10	ND	ND	ND	0.50	0	ND	0.56	104
196C Gamma-BHC	ND	.10	ND	ND	ND	0.18	88	ND	0.56	102
196D Delta-BHC	ND	.10	ND	ND	ND	0.50	164	ND	0.56	98
214 Kepone	ND	1.0	ND	ND	ND	5.0	70	ND	5.6	70
230 Methoxychlor	ND	10	ND	ND	ND	5.0	112	ND	56	106
242 Methyl parathion	ND	1.0	ND	ND	ND	5.0	84	ND	5.6	88
283 Parathion	ND	1.0	ND	ND	ND	5.0	96	ND	5.6	96
295 Phorate	-	-	-	-	-	-	-	-	-	-
296 Famphur	ND	10	ND	ND	ND	50	22	ND	56	27
333 Tetraethylpyrophosphate	ND	5.2	ND	ND	ND	25	48	ND	28	47
351 Toxaphene	ND	2.5	ND	ND	ND	20	115	ND	22	82
361 2,4,5-T	ND	5.0	ND	ND	ND	10	158	ND	10	104
362 2,4,5-TP (Silvex)	ND	5.0	ND	ND	ND	10	143	ND	10	103

*See Grand Table for complete compound name and specifications.

**EPR/ETC established method detection limit.

***"ND" specifies no standard available.

****"MDL" = Method Detection Limit. Major spikes could not be detected at current method levels.

SEP 30, 1985

TABLE 1: QUANTITATIVE RESULTS and QUALITY ASSURANCE DATA
RCRA Appendix 8 - DAI HPLC Compounds - Analysis Data (QR25)

Chain of Custody Data Required for ETC Data Management Summary Reports										
J3284 CBS RECORDS		CBSCARGWM WB9		850723 1030						
ETC Sample No.	Company	Facility	Sample Point	Date	Time	Elapsed Hours				
Compound *	Results	QC Replicate	QC Blank and Spiked Blank			QC Matrix Spike	Unspiked Sample mg/l	Concen. Added mg/l	Concen. Added mg/l	% Recov
	Sample Concen. mg/l	MDL** mg/l	First mg/l	Second mg/l	Blank Data mg/l	Concen. Added mg/l	% Recov	Unspiked Sample mg/l	Concen. Added mg/l	% Recov
6 1-Acetyl-2-thiourea	ND	.5	ND	ND	ND	5.0	90	ND	5.0	0
8 Acrylamide	ND	.5	ND	ND	ND	5.0	96	ND	5.0	95
17 Amitrole	-	-	-	-	-	-	-	-	-	-
82 1-(o-Chlorophenyl)thiourea	ND	.5	ND	ND	ND	5.0	83	ND	5.0	87
145 Diethylstilbestesterol	-	-	-	-	-	-	-	-	-	-
176 Ethyl carbamate	-	-	-	-	-	-	-	-	-	-
179 Ethyleneimine	-	-	-	-	-	-	-	-	-	-
181 Ethylenethiourea	ND	.5	ND	ND	ND	5.0	98	ND	5.0	96
221 Maleic hydrazide	ND	.5	ND	ND	ND	5.0	120	ND	5.0	130
222 Malononitrile	-	-	-	-	-	-	-	-	-	-
229 Methomyl	ND	.5	ND	ND	ND	5.0	39	ND	5.0	37
231 2-Methylaziridine	-	-	-	-	-	-	-	-	-	-
253 Nicotinic acid	ND	.5	ND	ND	ND	5.0	84	ND	5.0	130
260 Nitroglycerin	-	-	-	-	-	-	-	-	-	-
268 N-Nitroso-N-Ethylurea	ND	.5	ND	ND	ND	5.0	-	ND	5.0	-
270 N-Nitroso-N-methylurea	ND	.5	ND	ND	ND	5.0	-	ND	5.0	-
292 N-Phenylthiourea	ND	.5	ND	ND	ND	5.0	70	ND	5.0	74
309 Reserpine	-	-	-	-	-	-	-	-	-	-
345 Thiourea	ND	.5	ND	ND	ND	5.0	84	ND	5.0	130

*See Ground Table for complete compound name and specifications.
**EPA/ETC established method detection limit.
***specifies no standard available.
****Indeterminate. Standard and/or spikes could not be detected at current method levels.

SEP 30, 1985

TABLE 1: QUANTITATIVE RESULTS and QUALITY ASSURANCE DATA

RCRA Appendix 8 – Extractable HPLC Compounds – Analysis Data (QR31)

Chain of Custody Data Required for ETC Data Management Summary Reports									
J3284 CBS RECORDS		CBSCARGWM WB9		850723 1030					
ETC Sample No.	Company	Facility	Sample Point	Date	Time	Elapsed			
Compound *	Results	QC Replicate	QC Blank and Spiked Blank	QC Matrix Spike					
3 Warfarin	ND	25	ND	ND	500	120	ND	500	64
15 Mitomycin C	-	-	-	-	-	-	-	-	-
26 Azaserine	-	-	-	-	-	-	-	-	-
35 Benzidine	ND	25	BMDL	ND	ND	500	140	ND	500
62 Chlorambucil	-	-	-	-	-	-	-	-	-
86 Citrus red No. 2	-	-	-	-	-	-	-	-	-
99 Daunomycin	-	-	-	-	-	-	-	-	-
119 3,3'-Dichlorobenzidine	ND	25	BMDL	BMDL	BMDL	500	120	BMDL	500
147 3,4-Dihydroxy-alpha-(methyl-	-	-	-	-	-	-	-	-	-
150 3,3'-Dimethoxybenzidine	ND	25	BMDL	ND	ND	500	100	BMDL	500
153 3,3'-Dimethylbenzidine	ND	25	ND	ND	BMDL	500	83	ND	500
249 1-Naphthyl-2-thiourea	ND	25	ND	ND	ND	500	48	ND	500
290A m-phenylenediamine	ND	25	ND	ND	ND	500	65	ND	500
290B o-phenylenediamine	ND	25	ND	ND	ND	500	46	ND	500
290C p-phenylenediamine	BMDL	25	BMDL	ND	ND	500	49	ND	500
320 Streptozotocin	-	-	-	-	-	-	-	-	-
343 Thioacetamide	-	-	-	-	-	-	-	-	-
344 Thiosemicarbazide	-	-	-	-	-	-	-	-	-
348A Toluene-2,4-Diamine	ND	25	ND	ND	ND	500	110	ND	500
369 Trypan blue	-	-	-	-	-	-	-	-	-

*See Grand Table for complete compound name and specifications.
**EPA/ETC established method detection limit.

"ND" = Indeterminate. Standard and/or spikes could not be detected at current method levels.

AUG 17, 1985

TABLE 1: QUANTITATIVE RESULTS and QUALITY ASSURANCE DATA

RCRA Appendix 8 - Aroclors -GC Analysis Data (QR30)

Chain of Custody Data Required for ETC Data Management Summary Reports										
J3284	CBS RECORDS	CBSCARGWM	WB9	850723 1030						
ETC Sample No.	Company	Facility	Sample Point	Date	Elapsed Time	Hours				

Compound *	Results		QC Replicate		QC Blank and Spiked Blank			QC Matrix Spike		
	Sample Concen. ug/l	MDL** ug/l	First ug/l	Second ug/l	Blank Data ug/l	Concen. Added ug/l	% Recov	Unspiked Sample ug/l	Concen. Added ug/l	% Recov
300A Aroclor 1242	ND	1	ND	ND	ND	0	-	ND	0	-
300B Aroclor 1254	ND	1	ND	ND	ND	0	-	ND	0	-
300C Aroclor 1260	ND	1	ND	ND	ND	2	122	31	2	0
300D Aroclor 1248	ND	1	ND	ND	ND	0	-	ND	0	-
300E Aroclor 1232	ND	1	ND	ND	ND	0	-	ND	0	-
300F Aroclor 1221	ND	1	ND	ND	ND	0	-	ND	0	-
300G Aroclor 1016	ND	1	ND	ND	ND	0	-	ND	0	-

* Spiked samples that contain compounds present at high levels do not provide valid spike recovery data.

**See Grund Table for complete compound name and specifications.
***MDL calculated for each sample matrix.

SEP 1, 1985

TABLE 1: QUANTITATIVE RESULTS and QUALITY ASSURANCE DATA

RCRA Appendix 8 - Metals and Cyanides - Analysis Data (QR 29)

Chain of Custody Data Required for ETC Data Management Summary Reports									
J3284	CBS RECORDS	CBSCARGMM	WB9	850723	1030				
ETC Sample No.	Company	Facility	Sample Point	Date	Elapsed Time	Hours			

Compound *	Results		QC Replicate		QC Blank and Spiked Blank			QC Matrix Spike		
	Sample Concen. ug/l	MDL ug/l	First	Second	Blank Data	Concen. Added	% Recov	Unspiked Sample	Concen. Added	% Recov
13A Aluminum	2400	100								
19A Antimony	ND	200								
21A Arsenic	ND	10								
27A Barium	49	3								
42A Beryllium	ND	3								
56A Cadmium	ND	10								
57A Calcium	19400	200								
84A Chromium	ND	30								
88A Copper	BMDL	10								
92A Cyanide, Total	<25	25								
210 Iron	2200	200								
216A Lead	ND	200								
225A Mercury	BMDL	1								
250A Nickel	BMDL	30								
280A Osmium	BMDL	600								
301A Potassium	1200	40								
314A Selenium	ND	5								
317A Silver	ND	10								
319A Sodium	2400	80								
321A Strontium	ND	100								
335A Thallium	ND	5								
372A Vanadium	20	20								
374A Zinc	50	10								

*See Grand Table for complete compound name and specifications.

September 1, 1985

TABLE 2: METHOD PERFORMANCE DATA
Surrogate Recovery Water- GC/MS Data (QR20)
Chain of Custody Data Required for ETC Data Management Summary Reports
J3284

ETC Sample No. Company Facility Sample Point Date Time Elapsed Hours

Compound	Amount Added ug	% Recovery	Control Limits *	
			Lower	Upper
VOLATILE FRACTION				
Toluene-D8	250	100	86	119
Bromofluorobenzene	250	99	85	121
1,2-Dichloroethane-D4	250	98	77	120
ACID FRACTION				
Phenol-D5	100	29	15	103
2-Fluorophenol	100	62	23	121
2,4,6-Tribromophenol	100	100	10	130
BASE/NEUTRAL FRACTION				
Nitrobenzene-D5	50	78	41	120
2-Fluorobiphenyl	50	112	44	119
Terphenyl-D14	50	100	33	128

* IFB EPA Control Limits.

ETC Sample #: J3284**METHOD 8610**


The presence of total aromatics by ultraviolet absorption is determined by Method 8610. The sample is prepared by Method 3560. Method 3560 may be summarized as follows: the sample is drawn into a 50 ml syringe and a precleaned reverse phase cartridge is fitted on the syringe. The sample is forced through the cartridge. This cartridge is repeated with an additional 50 ml of sample. The cartridge is washed with distilled water and then eluted first with hexane and then methanol. The sample extracts are screened from 220 nm to 310 nm using an ultraviolet/visible spectrophotometer. Samples that have an absorbance greater than 0.005 absorbance units (99% transmittance) in this region fail the screen.

Extract	Pass	Fail
Hexane		X
Methanol		X

METHOD 8630

The presence of certain Appendix VIII compounds is determined by Method 8630. Aqueous samples are prepared using Method 3560 using diethyl ether as the eluting solvent. Method 3560 may be summarized as follows: The sample is drawn into a 50 ml syringe and a precleaned reverse phase cartridge is fitted on the syringe. The sample is forced through the cartridge. This cartridge is repeated with an additional 50 ml of sample. The cartridge is washed with distilled water and then eluted with diethyl ether. Ether extracts are dried over sodium sulfate and then reacted with lithium aluminum hydride. Compounds are then hydrolyzed to their parent alcohols or amines by adding a small amount of water and then the dinitrophenyl esters are formed by adding dinitrobenzoyl chloride. Excess benzoyl chloride is removed by adding water and washing with base. The ether extract is removed and screened from 220 nm to 700 nm using an ultraviolet/visible spectrophotometer. Samples that have an absorbance greater than 0.005 absorbance units (99% transmittance) in this region fail the screen.

Extract	Pass	Fail
Diethyl Ether		X



HSC 9/09/85 - 10/1/85




FIGURE 2. Mass spectrum of the 100% relative abundance base peak from Figure 1.

mass spectrum showing relative abundance (100%) for various ions (Table I).

m/z	Relative abundance of peak	Mass number	Relative abundance of peak	Mass number
37	100.00	95	52.83	53.93
59	30.60	93	58.93	59.93
63	10.00	97	72.25	73.25
83	10.00	99	74.25	75.25
95	100.00	101	82.19	83.19
103	10.00	103	6.67	6.67
125	10.00	125	81.53	82.53
147	10.00	147	6.66	6.66
151	10.00	151	6.66	6.66

ANALYST: RICHARD J. LOCHER
INJECTOR LINE: 200[°]C, 0.00 ml/min

DETECTOR: 200[°]C, 0.00 ml/min
SAMPLER: 200[°]C, 0.00 ml/min

RICHARD J. LOCHER
RECEIVED 1/29/85
REVIEWED 2/25/85
APPROVED 3/09/85
RICHARD J. LOCHER
RECEIVED 1/29/85
REVIEWED 2/25/85
APPROVED 3/09/85

FILE NUMBER APR/VAA, REARER, A P-RFR SUB FTF FAIRBANKS STA: Clean 110
BPK AB 4221 8.18 min.




TABLE 2: METHOD PERFORMANCE DATA (QR21)

GC/MS Tuning Data - Bromofluorobenzene (BFB) for Volatiles Analysis

m/z	Ion Abundance Criteria	% Relative Abundance		
		Base Peak	Appropriate Peak	Status
50	15-40% of mass 95	22.44	22.44	Ok
75	30-60% of mass 95	50.49	50.49	Ok
95	Base peak, 100% relative abundance	100.00	100.00	Ok
96	5-9% of mass 95	8.67	8.67	Ok
173	Less than 1% of mass 95	0.00	0.00	Ok
174	Greater than 50% of mass 95	89.17	89.17	Ok
175	5-9% of mass 174	7.20	8.08	Ok
176	95-101% of mass 174	85.69	96.09	Ok
177	5-9% of mass 176	6.68	7.80	Ok

Injection Date: 08/01/85

Injection Time: 23:55

Run No: >D0076

Spectrum No: 119

Analyst:

D. Schaefer

Processor:

B. O'Farrell

QC Batch:

QV 3555

Samples:

J3098, J3702, J3284, J2883

J2258, J2257, J3097, J3100, J4040

JB

File >E0712 AP7/8 ON E, 850809 DFTPP, 50NG, ETC CALIB CP Scan 370
Bpk Ab 34088 9.06 min.




TABLE 2: METHOD PERFORMANCE DATA (QR23)

GC/MS Tuning Data - Decafluorotriphenylphosphine (DFTPP) for ~~Base/Neutral~~
Analysis

m/z	Ion Abundance Criteria	% Relative Abundance		
		Base Peak	Appropriate Peak	Status
51	30-60% of mass 198	55.03	55.03	Ok
68	Less than 2% of mass 69	.92	1.42	Ok
69	(reference only)	64.49	64.49	Ok
70	Less than 2% of mass 69	.43	.66	Ok
127	40-60% of mass 198	53.63	53.63	Ok
197	Less than 1% of mass 198	0.00	0.00	Ok
198	Base peak, 100% relative abundance	100.00	100.00	Ok
199	5-9% of mass 198	6.41	6.41	Ok
275	10-30% of mass 198	20.02	20.02	Ok
365	Greater than 1% of mass 198	1.87	1.87	Ok
441	Less than mass 443	6.64	73.09	Ok
442	Greater than 40% of mass 198	48.44	48.44	Ok
443	17-23% of mass 442	9.09	18.77	Ok

Injection Date: 08/10/85

Injection Time: 13:24

Run No: >E0712

Spectrum No: 370

Analyst:

Processor:

QC Batch:

Samples:

Thomas M. Power
QC3418
J2057, J2258, J3285, J3290
J22956, J3702, J3284

WT
8/13/85

File >E0728 AP7/8 ON E, 850909 DFTPP, SONG, ETC CALIB CPD

Scan 366
9.06 min.

Bpk Ab 33856




TABLE 2: METHOD PERFORMANCE DATA (QR23)

GC/MS Tuning Data - Decafluorotriphenylphosphine (DFTPP) for Base/Neutral Analysis

M/z	Ion Abundance Criteria	% Relative Abundance		
		Base Peak	Appropriate Peak	Status
51	30-60% of mass 198	43.72	43.72	Ok
68	Less than 2% of mass 69	.83	1.67	Ok
69	(reference only)	49.83	49.83	Ok
70	Less than 2% of mass 69	.33	.66	Ok
127	40-60% of mass 198	43.93	43.93	Ok
197	Less than 1% of mass 198	.67	.67	Ok
198	Base peak, 100% relative abundance	100.00	100.00	Ok
199	5-9% of mass 198	6.33	6.33	Ok
275	10-30% of mass 198	22.41	22.41	Ok
365	Greater than 1% of mass 198	2.10	2.10	Ok
441	Less than mass 443	9.39	71.16	Ok
442	Greater than 40% of mass 198	67.39	67.39	Ok
443	17-23% of mass 442	13.19	19.57	Ok

Injection Date: 08/11/85

Analyst:

Injection Time: 05:10

Processor:

Run No.: >E0728

QC Batch:

Spectrum No.: 366

Samples: J2257, 58, J3285 to 90,

J2256, J3702, J3284




TABLE 2: METHOD PERFORMANCE DATA (QR23)

GC/MS Tuning Data - Decafluorotriphenylphosphine (DFTPP) for Base/Neutral Analysis

m/z	Ion Abundance Criteria	% Relative Abundance			Status
		Base Peak	Appropriate Peak		
51	30-60% of mass 198	52.60	52.60	Ok	
68	Less than 2% of mass 69	.81	1.08	Ok	
69	(reference only)	75.19	75.19	Ok	
70	Less than 2% of mass 69	0.00	0.00	Ok	
127	40-60% of mass 198	52.48	52.48	Ok	
197	Less than 1% of mass 198	0.00	0.00	Ok	
198	Base peak, 100% relative abundance	100.00	100.00	Ok	
199	5-9% of mass 198	6.79	6.79	Ok	
275	10-30% of mass 198	22.46	22.46	Ok	
365	Greater than 1% of mass 198	2.27	2.27	Ok	
441	Less than mass 443	14.80	93.62	Ok	
442	Greater than 40% of mass 198	88.19	88.19	Ok	
443	17-23% of mass 442	15.81	17.93	Ok	

Injection Date: 08/11/85

Injection Time: 12:53

Run No: >J4687

Spectrun No: 2912

Analyst:

Processor:

QC Batch:

Samples:

Thomas M. Price
Thomas M. Price
QC3418 Rev
J2252, J2258, J3285-J3290
J3702, J3284, J2951




TABLE 2: METHOD PERFORMANCE DATA (QR23)

ICP/MS Tuning Data - Decafluorotriphenylphosphine (DFTPP) for Base/Neutral Analysis Derivatized Formaldehyde

m/z	Ion Abundance Criteria	% Relative Abundance		
		Base Peak	Appropriate Peak	Status
51	30-60% of mass 198	42.53	42.53	OK
69	Less than 2% of mass 69 (reference only)	0.00	0.00	OK
127	40-60% of mass 198	43.89	43.89	OK
197	Less than 1% of mass 198	0.00	0.00	OK
198	Base peak, 100% relative abundance	100.00	100.00	OK
199	5-9% of mass 198	6.28	6.28	OK
226	10-30% of mass 198	21.27	21.27	OK
365	Greater than 1% of mass 198	1.92	1.92	OK
441	Less than mass 443	8.65	75.35	OK
442	Greater than 40% of mass 198	60.55	60.55	X 8124
443	17-23% of mass 442	11.48	18.95	OK

Injection Date: 08/18/85

Injection Time: 14:20

Run No: >18928

Spectrum No: 1027

Analyst:

Processor:

QC Batch:

Samples:

W. Hilti
Scalpinski

QB3429

J0545, J3284, J4933-J4935

Methodology For Analysis of Appendix VIII Parameters

The methods employed in the analysis of specific compounds in Appendix VIII are based on the second edition (July 1982) of EPA Manual SW-846 "Test Methods for Evaluating Solid Waste". Since the manual is deficient in details on how the complete analysis of Appendix VIII compounds can be performed, we devised an analytical scheme that divides the compounds into the twelve categories as stated in the introduction. Within the nine categories of compounds that can be specifically analysed, not all the reference standards are readily available. We footnoted those unavailable compounds in our result tables. In the GC/MS analyses we were able to search for the specific compounds that did not have corresponding reference standards by using standard spectra from the literature. In other cases, GC and HPLC, where retention time characteristics are the only means to identify the compounds, we were unable to draw any conclusions on those compounds that corresponding reference standards are not available at the time of analysis.

In each category rigid compliance with the instrument parameters and performance criteria of established EPA methods, such as those in methods 624 and 625, was achieved before any sample analysis was initiated. (Similar criteria do not exist in SW846.) For GC/MS analysis where standards were available, identification was performed using relative retention times, the relative abundance of three characteristic ions and the abundance ratios. The entire mass spectrum was reviewed to confirm each identification. Quantitative analysis was performed using an internal standard with a single characteristic ion. When compounds without corresponding reference standard were identified, they were quantified assuming the same response factor as the internal standard.

Volatile Compounds by Purge and Trap GC/MS

For the analysis of purge and trap Volatile compounds, Methods 8240 and 5030 were used. The analysis can be summarized as follows: Helium is bubbled through a 5 ml water sample contained in a specially designed purging chamber at ambient temperature. The purgeable volatile organic compounds are transferred from the aqueous phase to the vapor phase. The vapor is swept through a sorbent column where the organic components are trapped. After the purge cycle is complete, the sorbent column is heated and backflushed with helium to desorb the organic purgeables onto a gas chromatographic column. The gas chromatograph is temperature programmed to separate the purgeable mixture. The separated purgeable components are then identified and quantitated using a computerized mass spectrometer.

Water Soluble Compounds by Direct Aqueous Injection GC/MS

For the analysis of water soluble compounds, 5 ul of aqueous sample was injected directly into the GC/MS system. The chromatographic column employed in the procedure is the same column used for purge and trap analysis. MS scanning was begun prior to sample injection to capture mass intensity data for early eluting compounds. The GC oven temperature program used was that specified in procedure 8240.

Extractable Acid, Base/Neutral and Pesticide Compounds by GC/MS

For the analysis of the Acid, Base/Neutral and Pesticide compounds in water, Methods 3510 and 8270 were used. The analysis can be summarized as follows: a measured volume of sample, approximately 1 liter, is extracted with an aliquot of methylene chloride without pH adjustment and then the sample is adjusted to a pH greater than 11 and extracted with another aliquot of methylene chloride. These two aliquots were combined. The pH of the sample is then adjusted to a value less than 2 and extracted with another aliquot of methylene chloride. A separatory funnel or continuous extractor is used to perform the extractions. The two extracts are dried and concentrated to a 1 ml final volume. The extracts are then combined just prior to injection into a GC/MS instrument.

GC Analysis of Herbicides and Pesticides

The methods employed in the analysis of your sample for herbicides and pesticides are established EPA methods taken from the "Manual of Analytical Methods for the Analysis of Pesticides in Humans and Environmental Samples," June, 1980 and methods 8080, 8140 and 8150.

The herbicide method can be summarized as follows: A measured volume of sample, approximately 500-1000 ml, to which sodium sulfate has been added, is acidified and extracted with methylene chloride. The methylene chloride extract is evaporated to dryness, and the residue is derivatized with diazomethane and injected into a gas chromatograph equipped with a ^{63}Ni electron capture detector.

The pesticide method can be summarized as follows: A measured volume of sample, approximately 500ml, is extracted with methylene chloride. The extract is dried and concentrated to a final volume of 1ml and injected into a gas chromatograph equipped with a ^{63}Ni electron capture detector and Flame Photometric detector in phosphorus mode.

Analysis of Metals

The determination of metals in aqueous samples is performed according to the methods published by EPA in "Methods for Chemical Analysis of Water and Wastes," EPA-600/4-79-020, March, 1979, and the inorganic methods in SW846. Arsenic, selenium and thallium are determined by furnace AA; silver, aluminum, barium, beryllium, boron, cadmium, calcium, chromium, copper, cobalt, iron, magnesium, manganese, molybdenum, nickel, lead, sodium, antimony, tin, titanium, vanadium, and zinc are determined by ICP emission. The determination of mercury is performed by cold vapor AA.

Analysis by HPLC

The analysis of thermally unstable and polar compounds are based on the HPLC methods 8320 and 8330. The compounds analyzed fall into two catagories: direct aqueous injection and Base/Neutral extractables. Twenty μl of the sample or extract is injected into an HPLC equipped with a reverse phase column. Gradient elution and UV detector at 210 and 250 nm are used.

Conventional

Total cyanide analysis is performed using Method 9010.

Summary of Quality Assurance/Quality Control Procedures (QA/QC)

ETC bases its quality assurance protocols on the following government guidelines:

- . "Handbook for Analytical Quality Control in Water and Wastewater Laboratories", EPA-600/4-79-019, March 1979;
- . National Enforcement Investigation Center Policies, and Procedures manual; EPA-330/9/79/001-R, October 1979;
- . the recommended guidelines for EPA Methods 624 and 625. (Federal Register, December 3, 1979; updated on October 26, 1984);
- . "Manual of Analytical Methods for the Analysis of Pesticides in Humans and Environmental Samples," EPA 600/8-80-038, June 1980; and
- . "Determination of 2,3,7,8-TCDD in Soil and Sediment" EPA, Region VII, Kansas City, September 1983.
- . Organic Analysis: Multi-media, Multi Concentration-IFB WA84-A267
- . Dioxin Analysis: Soil/Sediment Matrix; Multi-Concentration; Selected Ion Monitoring with Jar Extraction Procedure-IFB WA84-A002

However, we have modified our protocols to provide a higher level of QA/QC than the guidelines require. For example, we analyze a higher than required number of quality control samples and we pay especially careful attention to the certification of the "reference standard" compounds we use in analysis. Below are listed the key QA/QC elements for the methods we used.

Analysis of Volatile Organic Compounds by Gas Chromatography/Mass Spectrometry

- Each batch of 13 samples consists of 9 customer samples (at a maximum), one blank sample, one spiked blank, one spiked sample and one replicate sample. This amounts to a 30% quality control factor.
- Three surrogate compounds are added to each sample in the batch of 13.
- A blind quality control sample is introduced to the laboratory for analysis on a weekly basis.
- Each GC/MS is checked and retuned, if necessary, at the beginning of each day to ensure that its performance on bromofluorobenzene (BFB) meets the EPA criteria.
- A calibration curve for quantitation is prepared using a mixture of Volatile Organic Priority Pollutant "standards" at a minimum of 3 different concentrations and using a mixture of 3 internal standards at a constant concentration.
- The calibration curve is verified with a mixture of priority pollutant standards every day. If the response factors vary greater than 25%, the instrument must be recalibrated.

Analysis of Organic Compounds Extracted in Acid or Base/Neutral Solutions by Gas Chromatography/Mass Spectrometry

- Each batch of 20 samples consists of 16 customer samples (at a maximum), one blank sample, one spiked blank (for water matrices), one sample spiked with the priority

pollutant standard mixture and a duplicate customer sample. This amounts to a 20% quality control factor.

- Three surrogate compounds are added to each sample in the batch for Base/Neutral analysis.
- Three surrogate compounds are added to each sample in the batch for Acid analysis.
- A blind quality control sample is introduced to the laboratory for analysis on a weekly basis.
- Each GC/MS is checked and retuned, if necessary, at the beginning of each day to ensure that its performance on decafluorotriphenylphosphine (DFTPP) meets the EPA criteria.
- A calibration curve for quantitation is prepared using a mixture of standards composed of either the Organic Acid or Base/Neutral Extractable Compounds at a minimum of 3 concentrations and using five internal standards for quantitation.

Analysis of Metals

All Samples

- New standards are prepared for each batch of samples.
- Normal calibration is performed using a blank sample and four standards that have been through the sample preparation procedure. A regression analysis is used to construct the calibration curve.
- All EP Toxicity samples and all samples determined by furnace atomic absorption are calculated by the "method of additions".
- For each sample analysis that requires the use of the "method of additions" technique, a three point calibration is performed using U.S. EPA "Methods for Chemical Analysis of Water and Wastes, 1979". Results are obtained using linear regression analysis. Any regression with a coefficient of correlation below 0.990 is considered suspect, necessitating review of calibration data or sample re-analysis.
- In constructing the normal calibration curves the lowest concentration levels we use are values greater than or equal to 5 times the Instrumental Detection Limit (IDL).
- All calibration standards are analyzed in duplicate, at a minimum.
- Independent reference standards are used to check the accuracy of calibration standards.
- A check standard is analyzed every ten samples to validate the normal calibration curve.
- One customer sample out of every ten is analyzed in triplicate.

Homogeneous Samples (except for Mercury analysis)

Samples are analyzed in batches of 30 or less. For batches in which the sample matrices are homogeneous, the QC program is a minimum of 25% and consists of analyzing:

- 3 sets of triplicate analyses;
- 2 Replicate spikes;
- 1 independent reference standard;
- 4 Calibration standards (processed using the sample preparation method);
- 4 Calibration standards (without sample preparation); and

- 1 Reagent Blank.

Heterogeneous Samples (except for Mercury analysis)

Samples are analyzed in batches of 30 or less. For batches in which the sample matrices are heterogeneous, the QC program is a minimum of 35% and consists of analyzing:

- 3 sets of triplicate analyses;
- 2 Replicate spikes;
- 1 Replicate independent reference standards;
- 4 Calibration standards (processed using the sample preparation method);
- 1 Procedural Blank;
- 4 Calibration standards (without sample preparation); and
- 1 Reagent Blank.

Analysis of Mercury

To analyze samples for mercury we group them by matrix in batches of 30 or less. Our QC program is a minimum of 30% and consists of analyzing:

- each of the 30 customer samples in duplicate;
- 3 sets of triplicate analyses;
- 2 Replicate spikes;
- 2 Replicate independent reference standards;
- 10 Calibration standards (processed using the sample preparation method); and
- 2 Procedural Blanks.

Analysis of Pesticides, Herbicides and PCB's by Gas Chromatography

Pesticide, herbicide and PCB samples are grouped in batches of 16 customer samples or less according to the type of analysis to be performed. The QC program for each of these three types of analyses is a minimum of 20% and consists of analyzing:

- 1 procedural blank sample(a reagent blank is analyzed in the case of non-water matrices);
- 1 spiked blank sample (the spiked blank is eliminated in the case of non-water matrices);
- 1 replicate sample;
- 1 replicate spiked sample; and
- 1 known reference QC sample for at least each 100 samples analyzed.

The instrument is calibrated each run with three standards, and checked every 10 samples.

Analysis of Cyanides, Phenols, Fluoride, Chloride, Nitrate and Nitrite

- All parameters are analyzed using a Technicon Autoanalyzer II GT.

- 3 calibration standards are analyzed at the beginning and end of each batch.
- Each batch (up to 80 samples) consists of analyzing one blank, one spiked blank, one duplicate and spiked sample every 20 samples, and an EPA known reference sample.

Analysis of Total Organic Carbon (TOC)

TOC samples are analyzed on a daily basis with the number of samples analyzed per day dependent on the request for duplicate or quadruplicate analyses. The quality control program is designed to maintain the appropriate amount of QC and consists of the following elements:

- Daily instrument calibration
- One blank
- Standard recalibration every 10 samples
- Spiked samples at a low and high level
- Every sample is run in duplicate at a minimum

Analysis of Total Organic Halide (TOX)

- Blank reagent water for absolute carbon background must contain less than 5 ug/l of halide (as chloride).
- Using a trichlorophenol standard, the mean adsorption efficiency must be within +/- 15% of the standard value.
- Calibration standards are run every 10 samples.
- Every sample is run in duplicate at a minimum.

Analysis of 2,3,7,8-TCDD (Dioxin) by GC/MS (SIM)

- Each sample is dosed with a known quantity of $^{13}\text{C}_{12}$ -2,3,7,8-TCDD as internal standard and $^{37}\text{Cl}_4$ -TCDD as surrogate standard. The action limits for surrogate standard results is +/- 40% of the true value. Samples showing surrogate standard results outside of these limits are reextracted and reanalyzed.
- Two laboratory "method blanks" are run along with each set of 24 or fewer samples. The method blank is also dosed with the internal standard and surrogate standard.
- At least one per set of 24 samples is run in duplicate to determine intralaboratory precision.
- Qualitative Requirements. The following are met in order to confirm the presence of native 2,3,7,8-TCDD:
 - a. Isomer specificity must be demonstrated initially and verified once per 8-hour work shift. The verification consists of injecting a mixture containing TCDD isomers which elute close to 2,3,7,8-TCDD. The 2,3,7,8-TCDD must be separated from interfering isomers, with no more than 25% valley relative to the 2,3,7,8-TCDD peak.
 - b. The 320/322 ratio is within the range of 0.67 to 0.87.
 - c. Ions 320, 322, and 257 are all present and maximize together the signal to mean noise ratio must be 2.5 to 1 or better for all 3 ions.
 - d. The retention time is equal (within 3 seconds) the retention time for the isotopically labeled 2,3,7,8-TCDD.

e. At least one of the positives can be confirmed by obtaining partial scan spectra from mass 150 to mass 350. The partial scan guidelines are as follows:

- . the 320/324 ratio should be 1.58 +/- 0.16
- . the 257/259 ratio should be 1.03 +/- 0.10
- . the 194/196 ratio should be 1.54 +/- 0.15

- One sample is spiked with native 2,3,7,8-TCDD at a level of 1.0 PPB (for soil) for each set of 24 or fewer samples.
- In cases where no native 2,3,7,8-TCDD is detected, the actual detection limit is estimated and reported based on a signal to noise ratio of 2.5 to 1 at ions 320 and 322.
- For each sample, the internal standard is present with at least a 10 to 1 signal to noise ratio for both mass 332 and mass 334. Also, the internal standard 332/334 ratio must be within the range of 0.67 to 0.87.

Subcontractor QA/QC

Each subcontractor is required to maintain an appropriate level of quality control. To insure this, each subcontractor is required to submit to ETC the quality control data for all analyses it performs. This data is kept on file at ETC. In general, the amount of quality control required is one duplicate sample with one spiked sample for every ten analyses.

Chain-of-Custody

The chain-of-custody procedure is part of our quality assurance protocol. We believe our chain-of-custody record fully complies with the legal requirements of federal, state and local government agencies and of the courts of law. The record covers:

- labeling of sample bottles, packing the Sample Shuttle and transferring the Shuttle under seal to the custody of a shipper;
- outgoing shipping manifests;
- the chain-of-custody form completed by the person(s) breaking the Shuttle seal, taking the sample, resealing the Shuttle and transferring custody to a shipper;
- incoming shipping manifests;
- breaking the Shuttle's reseal;
- storing each labeled sample bottle in a secured area;
- disposition of each sample to an analyst or technician; and
- the use of the sample in each bottle in a testing procedure appropriate to the intended purpose of the sample.

The records show for each link in this process:


- the person with custody; and
- the time and date each person accepted or relinquished custody.

Appendix A
Mass Spectral Data
for
Quantitated Compounds

- 1) A total ion chromatogram for each sample analyzed by a GC/MS instrument.
- 2) A Quant report used by the analyst to determine qualitative and quantitative results of the compounds present.
- 3) A mass spectrum and a reference spectrum for each priority pollutant compound detected in the sample.

TOTAL ION CHROMATOGRAM

File >D1405 27.0-160.0 amu. 850927,D,HEATED P&T J3284W,0V3774,L,1,1



Date File: /D1405:1:06
Name: 850927,D,HEATED P&T
Misc: J3284W,0V3774,L,1,1,

Id File: DAP8DI:US
Title: DH4 HEATED 8 IC FILE
Last Calibration: 850928 20:56

Operator ID: JDE275
Quant Time: 850928 21:31
Injected at: 850927 02:12


QUANT REPORT

Operator ID: JQ6275
Output File: ^D1405::A0
Data File: >D1405::U6
Name: d50927.D,HEATED Part
Misc: J3284V,QU3774,L,1,1,

Quant Rev: 4 Quant Time: 850928 21:31
 Injected at: 850927 02:12
 Dilution Factor: 1.00

ID File: DAP8D1::US
Title: DAI APENDIX B ID FILE
Last Calibration: 890928 20:50

Category P.T. Search Area Date Units



QUANT REPORT

Operator ID: DE-3476
 Output File: ^D00081::AD
 Data File: >D00081::U3
 Name: AP8ZUDA, 850801, D
 Model: 032840, QU3555,1,5,6,

Quant Rev: 4 Quant Time: 850801 12:35
 Injected at: 850802 03:17
 Dilution Factor: 1.00

ID File: DAP8PT::US

Title: IDELLE, PURGEABLE PRIORITY POLLUTANTS, I
 Last Calibration: 850801 23:52

	Compound	R.T.	Scan#	Area	Conc	Units	Q
10	*2-Bromo-1-chloropropane	20.04	486	6.0268	200.00	NG	96
211	Ethylbenzene	22.78	684	5.118	4.45	NG 10	85
212	Methylene chloride	24.48	162	44.07	34.01	NG-21=13	95
260	Toluene	25.17	612	1925	3.74	NG	90
281	1,1,1-Trichloroethane	14.95	354	1495	6.82	NG 1	88
331	1,2-Dichloroethane-04	(SURR) 13.63	320	25659	145.61	NG	86
341	Toluene-08	(SURR) 24.98	612	113589	250.73	NG	92
351	p-Bromofluorobenzene	(SURR) 30.58	256	53153	247.06	NG	99
391	1,2-Dibromoethane	20.04	486	2203	11.42	NG	26
411	1,4-Dichloro-2-butene	24.01	582	3118	13.38	NG 0	100
481	Methyl Ethyl Ketone	14.05	333	569	18.56	NG	94
511	*1,4-Dichlorobutane	24.01	582	62260	200.00	NG	97


* Compound is ISFD

RTL # 6

Instrument ID: 850818 21:35
Runstart Time: 850818 21:35
Operator ID: TR9113

File #: 10E11 E118 10E810117ED GRIMMINGS
ID E11A: PDRM: 115
Title: 10E11 E118 10E810117ED GRIMMINGS

Name: 850818 #1 GAPPZDFP
MS#: 132984B DB34791153115
Date: 219934:1116



GLJAN REPORT

Operator ID: 10913
Instrument ID: 10913:001
Date: 10/10/90
Name: 850518 #1 GAPER
Model: 10913, 002, 003, 004, 005

Op Date: 8/27/90


Time: 10:15

Site: DIFFERENT INSTRUMENTS

Last Calibration: 850610 18:39

Compound	RT (min)	Retention#	Area	Conc	Limits	RT
1) *ACETONE-DNP-H	24.32	1195	160120	2.50	11.00	
2) FORMAL DEHYD-DNP-H	24.52	41441	448	0.02	1.00	
* Compound is DNP						X 821

TOTAL INK CHROMATOGRAM
FILE: E:\5072\456459.D RUN: TIC



Data File: 8508221114

Name: APPENDIX E, 850819

Mass: 142.0417, 003416.1, 0011.1

Id File: EAPHE:US

Title: APPENDIX E TIC


Last Calibration: 850810 19:41

Operator ID: S13562
Run Time: 850811 23:45
Injected at: 850811 04:22

PML#42

TOTAL ION CHROMATOGRAM

File >E0727 45.0-450.0 amu. AP7/8 ON E, 850809 J3284C ,QC3418,L,9
TIC



Data File: >E0727::U4

Name: AP7/8 ON E, 850809

Misc: J3284C ,QC3418,L,990 ,1

BTL#32

Id File: EBNA8::US

Title: BN & ACID IDFILE FOR APPENDIX 8

Last Calibration: 850810 19:53

Operator ID: SJ3562

Quant Time: 850811 05:04

Injected at: 850811 04:22

QUANT REPORT

Operator ID: SJ3562
Output File: E0227::AQ
Data File: >E0227::U4
Name: AP7/8 ON E, 850809
Misc: 032840 ,QC3418,L,990 ,1

Quant Rev: 4 Quant Time: 850811 23:45
Injected at: 850811 04:22
Dilution Factor: 1.00

BTL#32


ID File: EAP8E::US
Title: APPENDIX 8 IDFILE
Last Calibration: 850810 19:40

Compound	R.T.	Scan#	Area	Conc	Units	q
15) *d4-1,4-Dichlorobenzene	8.95	302	146867	40.00	UG/ML	87
16) *d8-Naphthalene	12.36	495	349615	40.00	UG/ML	94
25) p-Chloroaniline	12.86	523	9920	5.26	UG/ML	5502648
29) *d10-Acenaphthalene	17.72	798	225194	40.00	UG/ML	97
46) *d10-Phenanthrene	22.20	1052	511089	40.00	UG/ML	100

* Compound is ISTD

TOTAL ION CHROMATOGRAM

File >E0745 45.0-450.0 amu. AP7/8 ON E, 850809 J3284C ,QC3418,L,9
TIC



Data File: >E0745::U4
Name: AP7/8 ON E, 850809
Misc: J3284C ,QC3418,L,990 ,1

BTL#50

Id File: EBNA8::US
Title: BN & ACID IDFILE FOR APPENDIX 8
Last Calibration: 850810 19:53

Operator ID: TM0576
Quant Time: 850812 01:21
Injected at: 850812 00:16

QUANT REPORT

Operator ID: TM0576
Output File: ^J4707::AQ
Data File: >J4707::U4
Name: AP8/DER ON J, 850811
Disc: J3284C ,QC3418,L,990 ,1

Quant Rev: 4 Quant Time: 850812 06:37
Injected at: 850812 05:54
Dilution Factor: 1.00
BTL#21

ID File: AP8DER::US
Title: DER AP8 IDFILe
Last Calibration: 850811 16:37

Compound	R.T.	Scan#	Area	Conc	Units	q
*d12-Chrysene	28.13	278	8.5060	46.00	UG/ML	100

* Compound is ISTD

QUANT REPORT

Operator ID: TM0576
Output File: EE0245::AQ
Data File: EE0245::04
Name: AP2/8 UN E, 850809
Misc: J3284C ,QC3418,L,990 ,1


Quant Rev: 4 Quant Time: 850813 06:56
Injected at: 850812 00:16
Dilution Factor: 1.00
BTL#50

ID File: EAPRE::US
Tit.: APFFENIX R TDE 11 F
Last Calibration: 850810 19:40

Compound	R.T.	Scan#	Area	Conc	Units	Q
1) *d4-1,4-Dichlorobenzene	9.03	298	189600	40.00	UG/ML	89
16) *d8-Naphthalene	12.44	492	427959	40.00	UG/ML	94
29) *d10-Acenaphthylene	12.79	295	268808	40.00	UG/ML	96
46) *d10-Phenanthrene	22.28	1049	614350	40.00	UG/ML	100

* Compound is DFTD

TOTAL ION CHROMATOGRAM



Date File: > J4707::U4
Name: AP8/DER ON J, 850911
Misc: J3284C ,QC3418,L,992 ,
S11 #21

Id File: AP8DER::US
Title: DER AP8 IDFILE
Last Calibration: 850811 16:37

Operator ID: TMU576
Quant Time: 850812 06:37
Injected at: 850812 05:54

QUANT REPORT

Operator ID: SJ3562
 Output File: ^E0727::AQ
 Data File: >E0727::U4
 Name: AP7/8 ON E, 850809
 Misc: J3284C ,QC3418,L,990 ,1

Quant Rev: 4 Quant Time: 850811 05:04
 Injected at: 850811 04:22
 Dilution Factor: 1.00
 BTL#32

ID File: EBNAB::US
 Title: BN & ACID IDFILE FOR APPENDIX B
 Last Calibration: 850810 19:53

	Compound	R.T.	Scan#	Area	Conc	Units	q
1)	*d4-1,4-Dichlorobenzene	8.95	302	146857	40.00	UG/ML	87
10)	2-Fluorophenol	(SURR)	7.08	196	12084	6.06	UG/ML
12)	Phenol-D5	(SURR)	8.95	302	660	.22	UG/ML
13)	*d8-Naphthalene		12.36	495	349615	40.00	UG/ML
20)	Naphthalene		12.86	523	5157	.66	UG/ML
27)	*d10-Acenaphthalene		17.72	798	225194	40.00	UG/ML
43)	*d10-Phenanthrene		22.20	1052	511089	40.00	UG/ML
49)	Di-n-butyl phthalate		24.48	1181	9436	.76	UG/ML
51)	*D12-Chrysene		30.37	1514	202101	40.00	UG/ML
56)	bis(2-Ethylhexyl)phthalate		31.01	1550	1632	.34	UG/ML
64)	Terphenyl-D14	(SURR)	27.27	1339	1287	.19	UG/ML

* Compound is ISTD

QUANT REPORT

Operator ID: TM0576
 Output File: ^E0745::AQ
 Data File: >E0745::U4
 Name: AP7/8 ON E, 850809
 Misc: J3284C ,QC3418,L,990 ,1

Quant Rev: 4 Quant Time: 850812 01:21
 Injected at: 850812 00:16
 Dilution Factor: 1.00

BTL#50

ID File: EBNAB::US
 Title: BN & ACID IDFILE FOR APPENDIX 8
 Last Calibration: 850810 19:53


	Compound	R.T.	Scan#	Area	Conc	Units	q
1)	*d4-1,4-Dichlorobenzene	9.01	298	189601	40.00	UG/ML	89
2)	N-Nitrosodimethylamine	4.28	37	823	.53	UG/ML	100
7)	Nitrobenzene-D5	(SURR)	10.51	383	137351	39.67	UG/ML
10)	2-Fluorophenol	(SURR)	6.39	150	144467	56.08	UG/ML
10)	2-Fluorophenol	(SURR)	6.45	153	5367	2.08	UG/ML
12)	Phenol-D5	(SURR)	8.43	265	108976	28.50	UG/ML
12)	Phenol-D5	(SURR)	9.01	298	936	.24	UG/ML
13)	*d8-Naphthalene	12.44	492	427939	40.00	UG/ML	94
14)	2-Fluorobiphenyl	(SURR)	15.83	684	199514	37.72	UG/ML
15)	N-Nitrosodi-n-propylamine	10.51	383	20045	8.67	UG/ML	61
27)	*d10-Acenaphthalene	17.79	795	268808	40.00	UG/ML	94
36)	Acenaphthene	17.88	800	781	.10	UG/ML	97
38)	Diethyl phthalate	19.53	893	2838	.38	UG/ML	85
39)	Fluorene	20.25	934	7932	.87	UG/ML	82
40)	2,4,6-Tribromophenol	(SURR)	20.25	934	149307	99.69	UG/ML
43)	*d10-Phenanthrene	22.28	1049	614350	40.00	UG/ML	96
49)	Di-n-butyl phthalate	24.56	1178	19632	1.32	UG/ML	87
51)	*D12-Chrysene	30.45	1511	218907	40.00	UG/ML	98
56)	bis(2-Ethylhexyl)phthalate	31.09	1547	2872	.55	UG/ML	99
64)	Terphenyl-D14	(SURR)	27.37	1337	362271	50.24	UG/ML
64)	Terphenyl-D14	(SURR)	27.43	1340	6246	.87	UG/ML

* Compound is ISTD

Appendix A1

Gas Chromatographic Spectral Data for Quantitated Compounds

- 1) A reconstructed gas chromatogram for each sample analyzed by a GC instrument.
- 2) A reconstructed gas chromatogram for the appropriate standard compounds analyzed with the same GC under the same operating conditions.



Sample: J3284

Injected at 23:40:08 ON AUG 5, 1985

Raw File: J3284T Proc File: J3284S Method: AP8FPD


Operator: AM

Case No: QC3413FPD


Instrument ID: #5


ul Injected: 4.5

Column ID: FPD




Sample: J3284 Injected at 22:08:53 ON AUG 6, 1985
Raw File: PA2880 Proc File: PA2880 Method: AP8PST
Operator: AM Case No: QC3413
Instrument ID: #5 μ l Injected: 4.5
Column ID: 1.5%SP2250/1.95%SP2401






Sample: AP8PST-C Injected at 8:27:28 ON AUG 7, 1985
Raw File: RA2695 Proc File: PA2695 , Method: AP8PST
Operator: AM Case No: QC3413
Instrument ID: #5 ul Injected: 4.5
Column ID: 1.5%SP2250/1.95%SP2401




Sample: APB190-C Injected at 14: 00: 49 ON AUG 7, 1985
Raw File: RA2703 Proc File: PA2703 , Method: AP8PST
Operator: AM Case No: QC3413
Instrument ID: #5 ul Injected: 4.5
Column ID: 1.5%SP2250/1.95%SP2401




188246

A/D Counts

15832



Sample: HERB-B Injected at 18:18:35 ON AUG 14, 1985
Raw File: RA2827 Proc File: PA2827 Method: HERB1



Sample: 245T-5.0 Injected at 18:28:37 ON AUG 14, 1985
Raw File: RA2828 Proc File: PA2828 Method: HERB1

Appendix A1

Gas Chromatographic Spectral Data for Quantitated Compounds

- 1) A reconstructed gas chromatogram for each sample analyzed by a GC instrument.
- 2) A reconstructed gas chromatogram for the appropriate standard compounds analyzed with the same GC under the same operating conditions.

5896382

A/D Counts

15818
0.000

4.817
4.477
4.027
3.577

4.729
4.017
3.637
3.277
2.893
2.545

5.578
4.729
4.017
3.637
3.277
2.893
2.545

6.482
6.177
7.725
8.000

9.758
9.758

11.42
12.00

14.00
16.00


18.00
20.00
22.00
24.00

Time (Minutes)

Sample: J3284

Injected at 2:42:03 ON AUG 2, 1985

Raw File: RJ1644 Proc File: PJ1644 Method: PCBB0



HP 5880A SAMPLER INJECTION @ 16:24 AUG 2, 1985
SAMPLE #: ID CODE :
15 1016-0.5

Ar 1221 -0.5ug/ml

15:58 8/2/85


3.22

0.35 0.94
1.11 1.46
2.03 2.47
2.64 3.08
3.22 3.64


1.58

4.04


RT: STOP RIN




Sample: AR01232-0.5 Injected at 16:24:49 ON AUG 2, 1985
Raw File: RN0177 Proc File: PN0177 Method: PCBO




Sample: AR1242-0.5 Injected at 23: 50: 04 ON AUG 1, 1985
Raw File: RJ1638 Proc File: PJ1638 Method: PCBB0




Sample: AR1248-0.5 Injected at 23: 21: 32 ON AUG 1, 1985
Raw File: RJ1637 Proc File: PJ1637 Method: PCBBO



Sample: AR1254-0.5 Injected at 22:53:03 ON AUG 1, 1985
Raw File: RJ1638 Proc File: PJ1638 Method: PCBB0





Sample: AR1260-0.5 Injected at 23: 50: 32 ON JUL 31, 1985
Raw File: RJ1608 Proc File: PJ1608 Method: PCBB0

Appendix A2
Chromatographic Spectral Data
for
Quantitated Compounds

- 1) A reconstructed chromatogram for each sample analyzed by an HPLC instrument.
- 2) A reconstructed chromatogram for the appropriate standard compounds analyzed with the same HPLC under the same operating conditions.

Rawdata copied from DPU memory to
date01/07/02

21/01/1900 19:13

Integration & online plot of new analysis


VIFL 23 07/01

SAPDAI

01/01/1900

20:24

J3284




AREA%

#	SI	TIME [min]	AREA [mAU]	TYPE	AREA%	HEIGHT [mAU]	WIDTH [min]	dTIME [min]	QUOTIENT [area]
1	A	21.145	42.117	BB	8.421	7.219	0.085	0.000	A/B 0.140
2	C	21.206	36.084	BV	13.878	3.840	0.235	0.002	C/D 1.024
3	A	23.073	458.04	BP	91.579	3.361	1.688	0.003	A/B 0.161
4	C	23.912	377.10	VB	69.969	2.002	2.337		
5	B	27.099	17.243	PB	0.499	1.054	0.365		
6	B	42.955	66.529	PV	1.927	1.264	0.726		
7	B	43.798	19.142	VB	0.554	0.816	0.405	0.004	B/C 0.233
8	B	58.435	210.28	ABB	6.091	3.416	0.809		
9	D	58.485	126.60	ABB	23.490	2.688	0.687		

TOTAL AREA FOR SIGNAL A = 500

TOTAL AREA FOR SIGNAL B = 3452

Integration & online plot of new analysis
UTAI 10 TNJ 1 RAPNAT
STB 11000-25



ARFA%

#	RT	TIMF min	ARFA min	TYPF min	ARFA%	HFTHT min	WTHTH min	HTMF min										
1	R	1.005	10.76R	-	RV	0.107	4.260	0.109	0.007	R/C	1.482	-	-	-	-	-	-	-

TELEGRAM REPORT

STD 11416 - 25

1	A	-24.100	210.4	270.4	252.1	260.1
2	A	10	10	10	10	10
3	RANGE					
4	RANGE					

5

6

7

8

9

10

11

12

13

14

15

16

17

18

19

20

21


22

23

24

25

PEAKWIDTH = 0.750



#	ST	TIME [min]	ARFA	TYPE	ARFA%	HFTANT [mAUs]	WIDTH [min]	RTMF [min]	QUINTENT [arcsec]	
1	A	21.153	15.439	RR	3.719	3.132	0.074	0.001	A/R	0.15?
2	A	21.244	7.921	RR	1.908	2.538	0.123	0.003	A/C	0.74?
3	D	22.234	63.777	UR	28.410	0.883	1.704			
4	A	23.091	2557.2	RU	59.190	20.094	1.497			
5	A	24.629	102.45	UR	3.573	2.410	0.709			
6	D	24.658	110.74	RR	2.768	4.077	0.391			
7	C	25.262	759.21	PU	46.680	1.148	3.995	0.010	A/R	0.369
8	A	25.269	31.931	UR	7.692	0.559	0.957	0.003	A/R	0.767
9	D	25.270	54.701	UR	6.470	0.877	0.979			
10	D	25.271	51.41	EP	5.759	2.987	0.965	0.002	R/C	1.849
11	D	25.272	14.10	RR	52.473	2.998	0.999			
12	D	25.273	14.10	RR	52.473	2.998	0.999			
13	D	25.274	14.10	RR	52.473	2.998	0.999			
14	D	25.275	14.10	RR	52.473	2.998	0.999			
15	D	25.276	14.10	RR	52.473	2.998	0.999			
16	D	25.277	14.10	RR	52.473	2.998	0.999			
17	D	25.278	14.10	RR	52.473	2.998	0.999			
18	D	25.279	14.10	RR	52.473	2.998	0.999			
19	D	25.280	14.10	RR	52.473	2.998	0.999			
20	D	25.281	14.10	RR	52.473	2.998	0.999			
21	D	25.282	14.10	RR	52.473	2.998	0.999			
22	D	25.283	14.10	RR	52.473	2.998	0.999			
23	D	25.284	14.10	RR	52.473	2.998	0.999			
24	D	25.285	14.10	RR	52.473	2.998	0.999			
25	D	25.286	14.10	RR	52.473	2.998	0.999			

Data copied from QPC memory to
date 01/01/1900

01/01/1900

08:41

Integration & online plot of new analysis


VIAL 00 INJ 1 SAPPDAI

01/01/1900

08:51

STD 11021-25 5.0

Signal	R	S	C	D	E	Annotation
Baseline	254.100	210.4	230.4	254.4	280.4	
Range	10	10	10	10	10	
Zero %	5	5	5	5	5	



PEAKWIDTH= 0.750

AREA%


#	SI	TIME [min]	AREA [mAU]	TYPE	AREA%	HEIGHT [mAU]	WIDTH [min]	dTIME [min]	QUOTIENT [area]
1	A	2.213	115.78	BB	6.438	9.203	0.184	0.000	A/B 0.349
2	A	3.447	248.08	BB	13.794	17.790	0.213	0.000	A/B 0.561
3	B	7.421	286.54	BV	4.159	10.451	0.428	0.002	B/C 4.970
4	A	9.387	213.99	VU	11.899	5.318	0.612	0.001	A/B 0.451
5	A	10.239	145.21	V8	8.074	3.622	0.617	0.000	A/B 0.579
6	C	21.195	26.243	BB	1.157	2.267	0.233		
7	A	23.068	471.77	BV	26.233	3.473	1.720	0.011	A/B 0.153
8	D	23.931	351.08	BB	24.376	1.708	2.600		
9	A	26.521	214.13	VU	11.906	8.426	0.399	0.001	A/B 0.374
10	A	28.374	65.072	VU	3.618	3.176	0.327	0.005	A/B 0.171

Integration & online plot of new analyses

SIGNAL A SIGNAL B SPECTRUM

STD 11417-05.0

Signal	A	B	C	D	E	Annotation
S. width	254,100	210,4	230,4	254,4	280,4	
Range	10	10	10	10	10	
Zero %	5	5	5	5	5	



AREA%

#	SI	TIME [min]	AREA [mAU]	TYPE	AREA%	HEIGHT [mAU]	WIDTH [min]	dTIME [min]	QUOTIENT [area]
1	B	21.197	15.105	BP	0.320	3.525	0.069		
2	A	21.197	445.32	BV	73.752	5.111	1.671	0.002	A/B 2.425
3	B	23.061	2501.7	VP	52.995	20.770	1.586		
4	D	23.942	517.52	BB	69.401	2.539	2.508		
5	B	27.228	20.768	PB	0.440	1.189	0.289		
6	B	31.654	272.57	VP	5.774	4.133	0.855	0.024	B/C 1.367
7	B	50.455	1359.5	BB	28.799	2.046	7.803		
8	C	56.368	4.661	BB	2.163	0.310	0.302		
9	B	58.315	367.39	ABB	7.783	6.005	0.800		
10	A	58.335	158.49	ABB	26.248	2.780	0.834		
11	D	58.360	222.82	ABB	29.881	4.004	0.767	0.005	D/E 3.197

TOTAL AREA FOR SIGNAL A = 604

TOTAL AREA FOR SIGNAL B = 4721

Rawdata copied from DDU memory to
Date 01/01/2000

01/01/2000

21:25


Integration & online plot of new analysis

VIAL 30 INJ 1 SAPDAI

01/01/1900

21:25

STD 11022 - 1.0



AREA%

#	SI	TIME [min]	AREA [mAU]	TYPE	AREA%	HEIGHT [mAU]	WIDTH [min]	dTIME [min]	QUOTIENT [area]
1	B	2.195	59.321	VB	1.582	4.571	0.178	0.001	B/C 1.231
2	A	3.450	52.376	BB	11.170	4.367	0.182	0.000	A/B 0.557
3	B	8.343	92.138	BB	2.457	2.523	0.568	0.000	B/C 0.728
4	C	21.188	18.311	BB	2.125	1.703	0.233	0.006	C/D 2.202
5	A	23.044	374.50	BV	79.869	3.017	1.603	0.017	A/B 0.141
6	A	26.521	42.018	VU	8.961	2.045	0.459	0.001	A/B 0.338
7	B	27.298	51.649	VU	1.377	1.525	0.631		
8	B	28.379	64.146	VU	1.711	1.666	0.540	0.005	B/C 1.400
9	B	29.084	127.71	VP	3.406	5.069	0.378	0.005	B/C 3.396
10	B	44.676	242.87	BB	6.477	0.829	3.557		
11	C	48.400	394.95	BB	45.843	0.528	8.891		

Data copied from CFS memory to
date 01/01/9000

01-01-1900 00:00

Integration & online plot of new analysis


VIAL 31 INJ 1 &APDRI

01/01/1900

00:00

STD11418-1.0

Signal	A	B	C	D	E	Annotation
3 wl,bw	254,100	210,4	230,4	254,4	280,4	
Range	10	10	10	10	10	
Zero %	5	5	5	5	5	



AREA%

#	SI	TIME [min]	AREA [mAU]	TYPE	AREA%	HEIGHT [mAU]	WIDTH [min]	dTIME [min]	QUOTIENT [area]
1	A	21.188	30.647	BV	10.904	2.639	0.230	0.004	A/C 2.660
2	A	23.012	250.42	VP	89.096	2.452	1.315		
3	B	23.051	2571.7	BP	84.913	21.005	1.544		
4	B	27.233	21.790	PB	0.719	1.212	0.335		
5	B	42.886	15.310	PB	0.506	0.755	0.348		
6	B	50.196	144.51	BP	4.771	1.824	1.100		
7	C	57.745	8.381	BP	42.106	0.361	0.387		
8	B	58.452	275.32	ABB	9.091	3.910	0.917		
9	D	58.489	131.71	ABB	96.252	2.835	0.700		

Integration & online plot of new analysis

VIAL 39 INT 1

DATAFILE

01:32:1303

30:12

Signal A 254,100 C 230,4

D 254,4


E 280,4

Annotation

CWL,BW 10,10 Range 10 Zero 5 STD 11023-0.5

dtIME 0.003 A/B 0.146

PEAKWIDTH= 0.750



AREA%

#	SI	TIME [min]	AREA [mAU]	TYPE	AREA%	HEIGHT [mAU]	WIDTH [min]	dTIME [min]	QUOTIENT [area]
1	B	3.367	46.216	BB	1.355	2.809	0.235	0.001	B/C 0.642
2	A	22.939	357.81	BU	100.00	2.610	1.586	0.003	A/B 0.146
3	D	23.951	184.93	BB	60.775	0.743	3.037		
4	B	26.333	44.397	UP	1.302	1.506	0.504		
5	C	26.353	19.132	BP	16.580	0.744	0.460		
6	B	27.161	13.372	PB	0.392	0.739	0.373		
7	B	49.844	618.16	BB	18.129	1.327	5.475		
8	C	57.781	24.264	BB	21.028	0.387	0.795		
9	B	58.478	223.64	ABB	6.735	3.156	0.971		
10	D	58.527	119.35	ABB	39.225	2.405	0.834		

Integration 3 online plot of new analysis
Run 42 INJ 1

STD 11419-0.5

Annotation

23:23

Signal	A	B	C	D	E
Chw1,bw	254,120	210,4	232,4	254,4	280,4
Range	10	10	10	13	13
Zero?	5	5	5	5	5

6

5

4


3

2

1

0

PEAKWIDTH= 0.750




AREA%

#	SI	TIME [min]	AREA [mAU]	TYPE	AREA%	HEIGHT [mAU]	WIDTH [min]	dTIME [min]	QUOTIENT [area]
1	A	22.933	344.91	BP	100.00	2.772	1.571	0.004	A/B 0.141
2	O	23.896	318.54	BB	72.052	1.459	2.985		
3	B	31.486	21.167	BB	0.774	0.793	0.477		
4	B	44.150	53.355	BB	1.951	0.524	1.373		
5	C	57.737	11.012	BP	100.00	0.380	0.483		
6	B	58.484	218.02	ABB	7.973	3.153	0.919		
7	D	58.536	123.56	ABB	27.948	2.340	0.732		

TOTAL AREA FOR SIGNAL A = 345
 TOTAL AREA FOR SIGNAL B = 2734
 TOTAL AREA FOR SIGNAL C = 11
 TOTAL AREA FOR SIGNAL D = 442

TOTAL AREA = 34530
TOTAL FID SIGNAL = 34530
DATA STREAM = 00000

#	RT (min)	AREA (COUNTS)	TYPE	AREA (COUNTS)	HEIGHT (mm)	WIDTH (min)	RT (min)	RT (min)
1	11.114	54.386	PB	0.157	3.329	0.132	0.860	1.070
2	11.427	5420.9	SP	14.243	1.365	0.593	0.303	0.378
3	11.469	172.38	CB	0.218	1.397	0.545	0.300	0.378
4	12.710	131.17	BD	0.134	0.597	0.501	0.300	0.378
5	13.328	475.37	BD	1.374	3.644	1.661	0.300	0.378
6	13.852	1245.6	PY	4.048	19.974	1.452	0.300	0.378
7	14.446	2567.8	YP	74.126	342.69	0.577	0.300	0.378
8	14.532	52.338	BD	0.863	1.739	0.513	0.300	0.378
9	14.703	362.54	BD	0.453	3.977	1.124	0.300	0.378
10	15.468	189.67	ACB	0.291	1.163	0.321	0.300	0.378
11	16.528	120.20	ABG	0.172	0.723	0.343	0.300	0.378



Ruth Sjöström

Qualitative Only

11/10/68

Nitroso methyl thiophene

3284
254.123
Range 75
Zero % 5

0 133.4
0 154.4
0 175.4
0 196.4


0 217.4
0 238.4
0 259.4

0 280.4
0 301.4
0 322.4

0 343.4
0 364.4
0 385.4

0 406.4

PEAKWIDTH= 3.753
Annotation



AREA%

#	SI	TIME [min]	AREA [mAU _s]	TYPE	AREA%	HEIGHT [mAU]	WIDTH [min]	dTIME [min]	QUOTIENT [area]
1	C	2.092	25.853	BV	0.577	3.924	0.097		
2	A	2.316	248.06	BV	8.605	31.499	0.106	0.001	A/B 0.191
3	B	2.575	36.587	VU	0.093	5.065	0.102		
4	A	2.732	23.812	VU	0.826	2.407	0.148	0.008	A/B 0.311
5	A	4.498	76.981	VU	2.670	18.528	0.064	0.000	A/B 0.201
6	B	19.928	363.10	BP	0.924	4.952	1.137		
7	A	21.750	257.51	BV	8.933	10.551	0.500	0.005	A/C 7.112
8	B	22.176	122.9.9	PV	3.129	30.569	0.600		
9	A	23.706	263.83	VU	9.152	5.335	0.824	0.007	A/B 0.177
10	A	24.219	467.40	UP	16.213	17.780	0.395	0.001	A/B 0.156
11	B	27.376	37.278	BP	6.695	2.648	0.273	0.004	B/C 0.834
12	A	28.779	542.11	PP	18.805	26.936	0.326	0.000	A/B 0.130
13	B	30.379	157.92	PV	0.402	7.221	0.386		
14	B	30.472	185.16	VU	6.471	8.910	0.346		
15	A	31.378	73.811	VP	2.560	2.196	0.436	0.004	A/B 0.431
16	A	32.503	487.97	PP	16.927	27.028	0.279	0.000	A/B 0.213
17	D	33.272	47.747	BP	4.524	2.611	0.416		

Integration & online plot of new analysis


VIAL 20 INJ 1 &APEXT

01/07/1900

23:31

STD 10897-1000

Signal	A	B	C	D	E	Annotation
S wl,bw	254,100	210,4	230,4	254,4	280,4	
Range	75	75	75	75	75	
Zero %	5	5	5	5	5	



AREA%

*	SI	TIME [min]	AREA [mAU]	TYPE	AREA%	HEIGHT [mAU]	WIDTH [min]	dTIME [min]	QUOTIENT [area]
1	B	2.170	12.363	VP	0.003	3.329	0.062		
2	A	3.480	80.592	BB	0.079	7.582	0.144	0.000	A/B 0.327
3	A	22.170	340.44	BV	0.332	30.395	0.187	0.001	A/B 0.231
4	A	22.708	476.54	VU	0.465	19.642	0.404	0.000	A/C 1.355
5	B	22.977	1485.5	VU	0.325	37.173	0.550		
6	C	23.714	295.72	VU	0.178	5.884	0.357		
7	A	27.215	768.88	VU	0.762	2.362	0.357		

dat 2011-07-02

INTEGRATION 3: CHROMATOGRAM ANALYSIS

FILE 21 INJ 1 ANALYSIS

STD 10898-1000

B 254,4

C 230,4

D 254,4

E 230,4

F 254,4


G 230,4

H 254,4

I 230,4

J 254,4

K 230,4




PEAKWIDTH= 3.750

AREA%							
#	SI	TIME [min]	AREA [mAU]	TYPE	AREA%	HEIGHT [mAU]	WIDTH [min]
1 A	2.152	80.786	BV	0.168	11.455	0.116	0.001 A/B
2 A	2.357	97.877	VV	0.204	9.157	0.152	0.011 A/B
3 A	2.550	114.33	VV	0.238	11.960	0.135	0.002 A/B
4 A	3.370	454.20	VB	0.947	85.550	0.080	0.000 A/B
5 D	21.844	161.24	BV	0.459	19.127	0.140	0.000 B/E
6 A	22.421	581.26	BV	1.212	13.091	0.655	0.008 A/C
7 B	22.756	1799.1	BV	0.409	28.462	0.307	0.366
8 D	23.426	78.958	VV	0.225	3.591	0.366	0.366
9 B	23.735	1725.8	VV	0.392	35.341	0.814	0.814
10 A	24.271	6192.0	VV	12.724	212.62	0.426	0.001 A/B
11 D	25.769	301.61	VV	0.858	8.736	0.457	0.457
12 A	25.792	213.90	VV	0.446	10.404	0.362	0.030 A/B
13 A	100.02						0.253 A/C


STD 10899-1080

Signal A 254,120 210,4 230,4 0 0
 S_w1,bw 15 75 75 5 5
 Range 15 75 75 5 5
 Terso % 5 5 5 5 5

Annotation



PEAKWIDTH = 0.753



AREA%


#	SI	TIME [min]	AREA [MAU]	TYPE	AREAT	HEIGHT [MAU]	WIDTH [min]	DTIME [min]	QUOTIENT [area]
1	B	2.163	9.432	BP	0.002	2.341	0.067		
2	A	3.433	15.113	BB	0.008	3.477	0.095	0.000	A/B 0.451
3	A	15.746	12243	BV	6.261	58.473	2.648	0.005	A/C 2.370
4	B	21.859	536.97	VV	0.088	20.685	0.435	0.019	B/C 0.312
5	A	21.900	1755.7	VV	0.898	57.132	0.491	0.004	A/D 0.521
6	C	23.375	63.666	VB	0.024	1.259	0.829		
7	A	23.553	252.79	VB	0.129	4.318	0.836		
8	B	23.591	1288.7	PP	0.212	23.702	0.891		
9	A	27.441	20.826	PV	0.011	2.744	0.215	0.004	A/E 0.373
10	A	28.757	57.372	CVV	0.029	2.230	0.429	0.001	A/B 0.362
11	A	29.410	23411	CVV	11.973	501.36	0.485	0.001	A/B 0.327
12	A	29.417	11111	CVV	1.987	0.419	0.007	0.007	A/B 0.312

Integration 3 - online plot of new analysis

10900 10900 10900

STD 10900-10900

Signal	A	B	C	D	E	Annotation
S, wl,Bw	254.100	210.4	233.4	254.4	200.4	
Range	75	75	75	75	75	
Zero %	5	5	5	5	5	



PEAKWIDTH= 3.750

AREA%

#	SI	TIME [min]	AREA [mAU]	TYPE	AREA%	HEIGHT [mAU]	WIDTH [min]	dTIME [min]	QUOTIENT [area]
1	A	2.425	15.957	BB	0.018	2.052	0.117	0.000	A/B 0.237
2	A	3.457	13.898	BB	0.016	3.159	0.115	0.001	A/B 0.493
3	A	21.841	219.15	BV	0.246	11.795	0.479	0.004	A/C 5.257
4	B	22.458	1026.4	BV	0.307	20.135	0.859		
5	A	23.724	225.94	VV	0.254	4.195	0.724	0.001	A/B 0.124
6	B	25.674	112.75	VV	0.034	3.895	0.447		
7	A	29.524	231.43	VV	0.260	3.438	0.892	0.004	A/B 0.395
8	D	29.790	1.937	BP	0.003	0.908	0.036		
9	D	35.468	1806.7	PV	2.923	19.286	1.561	0.001	D/E 1.333
10	A	36.174	41728	OVV	46.833	2263.6	0.281	0.001	A/S 0.338
11	C	36.243	72462	OVV	41.366	3599.3	0.336		
12	A	38.492	16595	VV	18.624	365.29	0.757	0.002	A/B 0.116
13	A	39.158	17744	VV	19.915	137.27	1.617	0.001	A/B 0.100
14	A	42.424	909.93	VV	1.021	30.865	0.412	0.007	A/B 0.217
15	A	44.131	2193.9	VV	2.462	104.90	0.318	0.001	A/B 0.222

Rawdata copied from CPU memory to
dat031:0720

01/08/1988

12:13

Integration & online plot of new analysis

VIAL 30 INJ 1


BAPEXT

01/08/1988

12:13


STD 10905-250

Signal	A	B	C	D	E	Annotation
S w1,bw	254,100	210,4	230,4	254,4	280,4	
Range	75	75	75	75	75	
Zero %	5	5	5	5	5	



AREA%

#	SI	TIME [min]	AREA [mAUa]	TYPE	AREA%	HEIGHT [mAU]	WIDTH [min]	dTIME [min]	QUOTIENT [area]
1	A	3.367	16.044	BP	0.040	4.066	0.091	0.001	A/B 0.489
2	A	21.826	59.989	BV	0.151	7.618	0.131	0.002	A/C 0.904
3	A	22.428	189.83	VU	0.477	4.979	0.570	0.014	A/D 2.301
4	B	22.803	1089.4	BV	0.629	20.557	0.966		
5	A	23.736	218.79	VU	0.550	4.105	0.752	0.001	A/B 0.143
6	A	24.608	8533.2	OVV	21.446	331.46	0.411	0.000	A/B 0.136
7	B	28.171	271.70	VU	0.157	3.432	1.319		
8	B	29.098	105.00	VU	0.061	2.651	0.613		
9	D	29.734	31.469	BP	0.158	1.055	0.492		
10	A	29.762	109.23	VU	0.275	1.449	0.979	0.007	A/C 0.316
11	B	29.805	176.20	VU	0.102	0.124	0.502		



#	SI	TIME [min]	AREA [mAU]	TYPE	AREA%	HEIGHT [mAU]	WIDTH [min]	dTIME [min]	QUOTIENT [area]
1	B	2.137	13.618	UV	0.024	4.247	0.053		
2	A	2.490	21.148	UV	0.170	2.462	0.126	0.000	A/B 0.458
3	A	3.220	11.741	UV	0.094	2.156	0.085	0.001	A/B 0.231
4	A	3.347	29.426	VB	0.237	8.222	0.060	0.000	A/B 0.528
5	A	21.818	63.253	BV	0.509	10.792	0.098	0.001	A/B 0.411
6	A	22.420	195.06	UV	1.569	4.603	0.616	0.006	A/C 2.217
7	B	22.788	934.58	UV	1.663	20.405	0.606		
8	A	23.741	213.16	UV	1.714	4.330	0.654	0.005	A/B 0.150
9	A	24.235	1407.8	VP	11.321	54.447	0.391	0.301	A/B 2.333
10	B	25.902	41.821	PV	0.074	1.728	0.410	0.018	B/C 0.525
11	B	27.684	32.386	VP	0.058	1.434	0.367		
12	A	28.414	4244.9	PV	34.136	113.71	0.516	0.000	A/B 2.173
13	A	29.342	1409.4	UV	11.334	31.508	0.596	0.001	A/B 0.152
14	B	32.584	2412.1	UV	4.291	38.763	0.807		
15	A	32.658	577.65	UV	4.645	8.739	0.956		

FILE NUMBER =

***** END OF REPORT *****

Rawdata copied from DPU memory to
dat311:D702


01/08/1990 12:40

Integration & online plot of new analysis

VIAL 32 INJ 1 &APEXT
STD 10907-250

01/08/1990 13:52

Signal	A	B	C	D	E	Annotation
S wl,bw	254,100	210,4	230,4	254,4	280,4	
Range	75	75	75	75	75	
Zero %	5	5	5	5	5	



AREA%

#	SI	TIME [min]	AREA [mAU _s]	TYPE	AREA%	HEIGHT [mAU]	WIDTH [min]	dTIME [min]	QUOTIENT [area]
1	A	3.370	14.356	BB	0.035	3.913	0.061	0.000	A/B 0.917
2	A	19.744	2280.6	BV	5.611	20.613	1.524	0.000	A/C 3.848
3	C	21.855	567.05	VP	1.823	19.546	0.420		

Integration & online plot of new analysis


VIAL 33 INJ 1 CAPERT

03/27/1995

1:18

STD 10908-25D

Signal	A	B	C	D	E	Annotation
S wl,bw	254,100	210,4	230,4	254,4	200,4	
Range	75	75	75	75	75	
Zero %	5	5	5	5	5	



AREA%

#	SI	TIME [min]	AREA [mAUa]	TYPE	AREA%	HEIGHT [mAU]	WIDTM [min]	dTIME [min]	QUOTIENT [area]
1	A	3.374	14.587	BB	0.091	3.549	0.101	0.000	A/B 0.430
2	A	21.748	56.609	BV	0.354	9.863	0.096	0.003	A/B 0.045
3	A	22.749	187.24	VV	1.171	3.735	0.679		
4	A	23.723	275.14	VP	1.721	5.110	0.715	0.003	A/B 0.150
5	A	29.760	34.240	PP	0.214	1.009	0.521	0.000	A/B 0.377
6	D	35.074	2.856	BP	0.025	2.320	0.135		
7	A	36.128	11285	OPV	70.597	688.23	0.265	0.000	A/B 0.337
8	A	38.475	861.63	VV	5.390	65.490	0.217	0.000	A/B 0.355
9	A	42.636	1230.4	VV	7.698	9.207	1.698	0.005	A/B 0.195
10	B	43.307	2169.5	VV	3.077	53.665	0.674		
11	A	43.898	738.35	VV	4.619	30.362	0.354	0.001	A/B 0.134
12	A	44.055	154.15	VV	1.027	4.771	0.577	0.000	A/B 0.177


VIAL 47 INJ 1
STD 10901-25

&APEXT

09/23/1985

13:43

Signal A	3	3	3	3	3
C w1.bw	154.100	212.4	230.4	154.4	230.4
Range	75	75	75	75	75
Zero %	5	5	5	5	5



AREA%

#	SI	TIME [min]	AREA [mAU _s]	TYPE	AREA%	HEIGHT [mAU]	WIDTH [min]	dTIME [min]	QUOTIENT [area]
1	A	3.342	14.609	BB	0.029	4.806	0.051	0.000	A/B 0.472
2	A	21.768	65.542	BV	0.132	10.263	0.106	0.002	A/C 0.729
3	A	22.251	188.31	VV	0.379	7.724	0.442	0.001	A/C 1.369
4	A	22.781	98.501	VV	0.198	3.967	0.414	0.005	A/B 0.065
5	C	23.552	260.88	VV	0.358	3.462	1.056		
6	A	23.704	265.66	VV	0.534	4.955	0.871	0.004	A/B 0.141
7	A	24.560	1474.5	VV	2.967	65.203	0.358	0.000	A/B 0.186
8	B	26.941	394.24	VV	0.298	7.167	0.725		
9	A	30.397	327.40	VV	0.659	4.079	1.015	0.003	A/B 0.191
10	D	30.918	48.764	VV	0.263	2.304	0.353		
11	A	31.535	477.22	VV	0.960	7.326	0.904	0.003	A/B 0.373
12	A	32.220	566.05	VV	1.139	6.752	1.103	0.000	A/B 0.631
13	C	32.518	1717.4	VV	2.354	14.730	1.467		
14	B	33.134	605.50	VV	0.457	15.980	0.519		
15	D	33.173	67.402	VV	0.363	1.127	0.997		
16	A	33.746	240.92	VV	0.485	5.368	0.748	0.019	A/B 0.333

Integration & baseline plot of new analysis
VIAL 43 INJ 1

STD 10902-25

DAPCXT

33/12/1995

14:32

Channel A

0 154.120 210.4 230.4 254.4 280.4

75

5

A/B

0.382

Channel B

0 0 0 0 0 0

75

5

A/C

1.536

Channel C

0 0 0 0 0 0

75

5

A/D

0.032

Channel D

0 0 0 0 0 0

75

5

A/E

0.143

Channel E

0 0 0 0 0 0

75

5

A/F

0.447

Channel F

0 0 0 0 0 0

75

5

A/G

0.210

Channel G

0 0 0 0 0 0

75

5

A/H

0.243

Channel H

0 0 0 0 0 0

75

5

A/I

0.341

Channel I

0 0 0 0 0 0

75

5

A/J

0.371

Channel J

0 0 0 0 0 0

75

5

A/K

0.341

Channel K

0 0 0 0 0 0

75

5

A/L

0.371

Channel L

0 0 0 0 0 0

75

5

A/M

0.371

Channel M

0 0 0 0 0 0

75

5

A/N

0.371

Channel N

0 0 0 0 0 0

75

5

A/O

0.371

Channel O

0 0 0 0 0 0

75

5

A/P

0.371

Channel P

0 0 0 0 0 0

75

5

A/Q

0.371

Channel Q

0 0 0 0 0 0

75

5

A/R

0.371

Channel R

0 0 0 0 0 0

75

5

A/S

0.371

Channel S

0 0 0 0 0 0

75

5

A/T

0.371

Channel T

0 0 0 0 0 0

75

5

A/U

0.371

Channel U

0 0 0 0 0 0

75

5

A/V

0.371

Channel V

0 0 0 0 0 0

75

5

A/W

0.371

Channel W

0 0 0 0 0 0

75

5

A/X

0.371

Channel X

0 0 0 0 0 0

75

5

A/Z

0.371

Channel Z

0 0 0 0 0 0

75

5

A/Y

0.371

Channel Y

0 0 0 0 0 0

75

5

A/P

0.371

Channel P

0 0 0 0 0 0

75

5

A/Q

0.371

Channel Q

0 0 0 0 0 0

75

5

A/R

0.371

Channel R

0 0 0 0 0 0

75

5

A/S

0.371

Channel S

0 0 0 0 0 0

75

5

A/T

0.371

Channel T

0 0 0 0 0 0

75

5

A/U

0.371

Channel U

0 0 0 0 0 0

75

5

A/V

0.371

Channel V

0 0 0 0 0 0

75

5

A/W

0.371

Channel W

0 0 0 0 0 0

75

5

A/X

0.371

Channel X

0 0 0 0 0 0

75

5

A/Z

0.371

Channel Z

0 0 0 0 0 0

75

5

A/Y

0.371

Channel Y

0 0 0 0 0 0

75

5

A/P

0.371

Channel P

0 0 0 0 0 0

75

5

A/Q

0.371

Channel Q

0 0 0 0 0 0

75

5

A/R

0.371

Channel R

0 0 0 0 0 0

75

5

A/S

0.371

Channel S

0 0 0 0 0 0

75

5

A/T

0.371

Channel T

0 0 0 0 0 0

75

5

A/U

0.371

Channel U

0 0 0 0 0 0

75

5

A/V

0.371

Channel V

0 0 0 0 0 0

75

5

A/W

0.371

Channel W

0 0 0 0 0 0

75

5

A/X

0.371

Channel X

0 0 0 0 0 0

75

5

A/Z

0.371

Channel Z

0 0 0 0 0 0

75

5

A/Y

0.371

Channel Y

0 0 0 0 0 0

75

5

A/P

0.371

Channel P

0 0 0 0 0 0

75

5

A/Q

0.371

Channel Q

0 0 0 0 0 0

75

5

A/R

0.371

Channel R

0 0 0 0 0 0

75

5

A/S

0.371

Channel S

0 0 0 0 0 0

75

5

A/T

0.371

Channel T

0 0 0 0 0 0

75

5

A/U

0.371

Channel U

0 0 0 0 0 0

75

5

A/V

0.371

Channel V

0 0 0 0 0 0

75

5

A/W

0.371

Channel W

0 0 0 0 0 0

***** END OF REPORT *****


Rawdata copied from DRU memory to
da2431:0701

09/20/1985

14:57

Integration & online plot of new analysis

VIAL 40 INJ 1 31/FEXT 09/20/1985 14:57
STD 10903-25



AREA%

#	SI	TIME [min]	AREA [mAU _s]	TYPE	AREA%	HEIGHT [mAU]	WIDTH [min]	dTIME [min]	QUOTIENT [area]
1	A	3.333	11.614	BB	0.065	3.870	0.075	0.000	A/B 0.525
2	A	21.995	709.03	BV	3.953	25.029	0.405	0.005	A/C 3.122
3	B	22.730	1224.0	BV	2.425	21.858	0.726		
4	A	23.701	276.94	VP	1.544	4.796	0.768	0.005	A/D 0.182
5	A	29.381	736.17	BV	4.439	52.972	0.238	0.000	A/B 0.301
6	B	30.379	31.122	VP	0.062	1.912	0.252	0.018	B/C 1.125
7	B	30.929	27.188	PV	0.054	1.587	0.285	0.001	B/C 1.057
8	C	31.584	36.313	VV	0.072	1.160	0.522		
9	B	32.224	55.205	VV	0.109	2.405	0.345		

INTEGRATION 3-DIMENSIONAL INTEGRATION ANALYSIS REPORT

23/25/1985

Annotation

Signal A 254.100 212.4 230.4 254.4 222.4


Signal B 254.100 212.4 230.4 254.4 222.4

Signal C 254.100 212.4 230.4 254.4 222.4

Signal D 254.100 212.4 230.4 254.4 222.4

Signal E 254.100 212.4 230.4 254.4 222.4

PEAKWIDTH= 0.752



AREA%

#	SI	TIME [min]	AREA [MAU]	TYPE	AREA%	HEIGHT [MAU]	WIDTH [min]	dTIME [min]	QUOTIENT [area]
1	A	3.343	12.606	BB	0.609	4.139	0.093	0.000	A/B 0.438
2	A	21.745	51.469	BV	2.487	8.714	0.998	0.000	A/B 0.344
3	C	22.698	97.083	UV	3.536	1.987	0.695		
4	A	22.728	188.43	UV	9.104	3.743	0.698		
5	A	23.713	272.52	UP	13.167	4.845	0.801	0.002	A/B 0.154
6	B	30.430	36.861	PP	0.543	1.014	0.666		
7	B	30.918	14.070	PB	0.207	1.159	0.292	0.005	B/C 0.735
8	B	32.233	44.419	BB	0.655	1.689	0.477		
9	A	36.025	1480.5	PP	71.531	107.57	0.231	0.000	A/B 0.432
10	A	38.376	64.199	PP	3.102	5.989	0.179	0.002	A/B 0.352


TOTAL AREA FOR SIGNAL A = 2970
 TOTAL AREA FOR SIGNAL B = 6785
 TOTAL AREA FOR SIGNAL C = 2746
 TOTAL AREA FOR SIGNAL D = 1063
 TOTAL AREA FOR SIGNAL E = 752

MUL FACTOR = 1

AREA% FOR SIGNAL A

Appendix A3


Total Aromatics and Derivatized Appendix VIII Compounds
by
UV/VIS



J3284 / QC3426
11/16/2010


Method J810 - QC3426

ABSORBANCE



J3284 / QC3426

Method 8810 - Q5212



J3284/QC3511


4041.W 8.30-QC3511

Appendix C1
GC/MS Subsidiary Data

TOTAL ION CHROMATOGRAM

File >D1391 27.9-160.0 amu. 850925,D,AP8HP&T 003774V ,QU3774,S,1,

TIC



Data File: >D1391

Name: 850925,D,AP8HP&T

Misc: QU3774V ,QU3774,S,1,1,

Id File: DAP8DI::US

Title: DAI APENDIX 8 ID FILE

Last Calibration: 850926 06:25

Operator ID: RR0835

Quant Time: 850926 06:26

Injected at: 850925 21:30

QUEEN'S REPORT

Operator ID: RR0835
Output File: ^D1391::U6
Data File: >D1391::U1
Name: 8611925,U,AP8HP&T
Misc: QC3774U ,QU3774,S,1,1,


Quant Rev: 4 Quant Time: 8511926 06:26
 Injected at: 8511925 21:30
Dilution Factor: 1.00

ID File: DAP8DI::US
Title: DAI APENDIX 8 ID FILE
Last Calibration: 850926 06:25

Compound	R.T.	Scan#	Area	Conc	Units	q
----------	------	-------	------	------	-------	---

TOTAL ION CHROMATOGRAM

File >D0072 40.0-265.0 amu. AP8/VOA, 850801, D QC3555V ,QV3555,L,5,5,BL
TIC



Data File: >D0072::U1

Name: AP8/VOA, 850801, D

Misc: QC3555V ,QV3555,L,5,5,BLANK

Id File: DAP8PT::US

Title: IDFILE, PURGEABLE PRIORITY POLLUTANTS, I

Last Calibration: 850801 23:52

Operator ID: DS3476

Quant Time: 850802 00:55

Injected at: 850801 18:02

QUANT REPORT

Operator ID: DS3476
 Output File: ^D0072::AQ
 Data File: >D0072::U1
 Name: AP8/VDA, 850801, D
 Misc: QC3555V ,QV3555,L,5,5,BLANK


Quant Rev: 4 Quant Time: 850802 00:55
 Injected at: 850801 18:02
 Dilution Factor: 1.00

ID File: DAP8PT::US
 Title: IDFILE, PURGEABLE PRIORITY POLLUTANTS, I
 Last Calibration: 850801 23:52

Compound	R.T.	Scan#	Area	Conc	Units	q
1) *2-Bromo-1-chloropropane	20.00	484	45093	200.00	NG	97
23) Methylene chloride	7.52	163	2490	21.50	NG	84
26) Toluene	25.13	616	2233	4.83	NG	96
28) 1,1,1-Trichloroethane	14.95	354	1970	10.09	NG	90
33) 1,2-Dichloroethane-D4	(SURR) 13.59	319	25067	265.40	NG	89
34) Toluene-D8	(SURR) 24.94	611	111847	275.20	NG	93
35) p-Bromofluorobenzene	(SURR) 30.46	753	53029	274.75	NG	97
39) 1,2-Dibromoethane	20.00	484	2025	46.82	NG ^b	78
41) 1,4-Dichloro-2-butene	23.93	585	2426	42.44	NG ^b	100
53) *1,4-Dichlorobutane	23.96	586	51700	200.00	NG	97

* Compound is ISTD

TOTAL ION CHROMATOGRAM

File >E0713 45.0-450.0 amu. AP7/8 ON E, 850809 QC3418C ,QC3418,L,1
TIC

Data File: >E0713::U4
Name: AP7/8 ON E, 850809
Misc: QC3418C ,QC3418,L,1000,1

BTL#18

Id File: EBNAB::US
Title: BN & ACID IDFILE FOR APPENDIX 8
Last Calibration: 850810 19:53

Operator ID: TM0576
Quant Time: 850810 21:27
Injected at: 850810 13:44

QUANT REPORT

Operator ID: TM0576
 Output File: ^E0713::AQ
 Data File: >E0713::U4
 Name: AP7/8 ON E, 850809
 Misc: QC3418C ,QC3418,L,1000,1

Quant Rev: 4 Quant Time: 850810 21:27
 Injected at: 850810 13:44
 Dilution Factor: 1.00
 BTL#18


ID File: EBNA8::US
 Title: BN & ACID IDFIL FOR APPENDIX 8
 Last Calibration: 850810 19:53

	Compound	R.T.	Scan#	Area	Conc	Units	q	
1)	*d4-1,4-Dichlorobenzene	8.94	304	160660	40.00	UG/ML	91	
10)	2-Fluorophenol	(SURR)	6.34	157	69531	31.85	UG/ML	100
12)	Phenol-D5	(SURR)	8.34	270	77512	23.93	UG/ML	90
12)	Phenol-D5	(SURR)	8.94	304	677	21 UG/ML	850823	2
13)	*d8-Naphthalene	12.35	497	361522	40.00	UG/ML	94	
27)	*d10-Acenaphthalene	17.70	800	218733	40.00	UG/ML	94	
43)	*d10-Phenanthrene	22.19	1054	387224	40.00	UG/ML	97	
49)	Di-n-butyl phthalate	24.47	1183	11015	1.17	UG/ML	96	
51)	*D12-Chrysene	30.34	1515	103348	40.00	UG/ML	97	
56)	bis(2-Ethylhexyl)phthalate	30.99	1552	945	38 UG/ML	850823	93	

* Compound is ISTD

TOTAL ION CHROMATOGRAM

File >E0731 45.0-458.0 amu. AP7/8 ON E, 850809 QC3418C ,QC3418,L,1
TIC



Data File: >E0731::U4
Name: AP7/8 ON E, 850809
Misc: QC3418C ,QC3418,L,1000,1

BTL#36

Id File: EBNAB::US
Title: BN & ACID IDFILE FOR APPENDIX 8
Last Calibration: 850810 19:53

Operator ID: SJ3562
Quant Time: 850811 07:44
Injected at: 850811 07:01

QUANT REPORT

Operator ID: SJ3562
 Output File: ^E0731::AQ
 Data File: >E0731::U4
 Name: AP7/8 ON E, 850809
 Misc: QC3418C ,QC3418,L,1000,1

Quant Rev: 4 Quant Time: 850811 07:44
 Injected at: 850811 07:01
 Dilution Factor: 1.00
 BTL#36


ID File: EBNAB::US
 Title: BN & ACID IDFILE FOR APPENDIX 8
 Last Calibration: 850810 19:53

	Compound	R.T.	Scan#	Area	Conc	Units	q
1)	*d4-1,4-Dichlorobenzene	8.95	305	192120	40.00	UG/ML	88
2)	N-Nitrosodimethylamine	4.29	42	848	.54	UG/ML	100
7)	Nitrobenzene-D5	(SURR) 10.45	390	149358	42.58	UG/ML	97
10)	2-Fluorophenol	(SURR) 6.37	159	171218	65.60	UG/ML	100
12)	Phenol-D5	(SURR) 8.37	272	116749	30.14	UG/ML	92
12)	Phenol-D5	(SURR) 8.95	305	884	.23	UG/ML	850821
13)	*d8-Naphthalene	12.38	499	434221	40.00	UG/ML	95
14)	2-Fluorobiphenyl	(SURR) 15.77	691	209838	39.10	UG/ML	99
15)	N-Nitrosodi-n-propylamine	10.45	390	21999	9.38	UG/ML	85082361
27)	*d10-Acenaphthalene	17.73	802	270380	40.00	UG/ML	94
39)	Fluorene	20.19	941	6927	.76	UG/ML	8508237
40)	2,4,6-Tribromophenol	(SURR) 20.19	941	125334	83.20	UG/ML	95
43)	*d10-Phenanthrene	22.20	1055	505143	40.00	UG/ML	97
49)	Di-n-butyl phthalate	24.50	1185	10114	.83	UG/ML	8508237
51)	*D12-Chrysene	30.37	1517	171335	40.00	UG/ML	97
64)	Terphenyl-D14	(SURR) 27.29	1343	305086	54.06	UG/ML	100
64)	Terphenyl-D14	(SURR) 27.34	1346	4829	.86	UG/ML	850823100

* Compound is ISTD

TOTAL ION CHROMATOGRAM

File >J4691 45.0-450.0 amu. AP8/DER ON J, 850911QC3418C ,QC3418,L,10
TIC



Data File: >J4691::U6

Name: AP8/DER ON J, 850911

Misc: QC3418C ,QC3418,L,1000,1

BTL# 5

Id File: AP8DER::US

Title: DER AP8 IDFILE

Last Calibration: 850811 16:37

Operator ID: TM0576

Quant Time: 850811 17:16

Injected at: 850811 16:27


QUANT REPORT

Operator ID: TM0576 Quant Rev: 4 Quant Time: 850811 17:16
Output File: ^J4691::AQ Injected at: 850811 16:27
Data File: >J4691::U6 Dilution Factor: 1.00
Name: AP8/DER DN J, 850911
Misc: QC3418C ,QC3418,L,1000,1 BTL# 5

ID File: AP80ER::US
Title: DER AP8 IDFILE
Last Calibration: 850811 16:37

Compound	R.T.	Scan#	Area	Conc	Units	q
1) *d12-Chrysene	28.13	1252	42234	40.00	UG/ML	100

* Compound is ISTD



Data File: 1R932:016
 Name: R60518 #1 GAPS/DER
 Msc: 0R3429R,0R3429,I,S,O,S
 RT: # 4

Id File: FORM:016
 Title: INFIL F DRIVITIZFD UNMPIINDS
 Last Calibration: 860610 18:39

Operator ID: TR9113
 Quant Time: 860818 20:07
 Injected At: 860818 19:29

QUANT REPORT

Operator ID: TR9113
Output File: ^18932::AQ
Data File: >18932::116
Name: R50518 #1 GAPB/DFR
Misc: QC3429B,QA3429,I,5,0.5

Quant Rev: 4 Quant Time: 850818 20:02
Injected at: 850818 19:29
Dilution Factor: 1.000
RIN #: 4

ID File: FORM::US
Title: IDEFILE FOR DERIVATIZED COMPOUNDS
Last Calibration: 850610 18:39

Compound	R.T.	Scan#	Area	Conc	Units	%
1) *ACETONE-DNPH	24.37	1198	192373	2.50	UG/ML	27
2) <u>FORMAL DEHYDRO</u> -DNPH	24.40	1142	4320	.08	UG/ML	8124

* Compound is ISTD

Appendix E

Chain-of Custody Forms

- 1) A field Chain-of-Custody form (CC1) is included for all samples shipped by ETC shuttle.
- 2) An in-house sample Chain-of Custody form is included for all samples not shipped by ETC shuttle.
- 3) Any additional Chain-of-Custody material provided by a client or by a client's sampling agent is also included.
- 4) A subcontractor's Chain-of-Custody form is included for any analytical work not performed within ETC's laboratory.
- 5) Analysis and Extraction Custody forms are included for the period the sample was in ETC's possession.

**ETC ENVIRONMENTAL
TESTING AND CERTIFICATION**

CHAIN OF CUSTODY FORM (CC1)

Shuttle A

Seal No. 37077 ETC Job# 53284

Date Sealed 7/12/85 By: Vlora

Company: CBS Records

Attn.: mark donas

Facility/Site: _____

Phone: (404) 838-2343

Address: 5152 Columbia Drive

Carrollton, GA 30117

SAMPLE IDENTIFICATION

Facility: 10

Sample Point: W-B-9 **Your Sample Point ID:** 072385 **Start Date:** 10/30/00
Source Code (from below) **Start Time:** 12:00 PM (clock) **Elapsed Month:**
 (left justify) (MM/DD/YY) (composite)

Source Codes: Well ..(W) Outfall.....(O) Bottom Sediment(B) Surface Impoundment....(I) Leachate Collection Sys....(C) Other(X)
Soil ...(S) River/Stream... (R) Generation Point(G) Treatment Facility(T) Lake/Ocean(L) Specify _____

SHUTTLE CONTENTS

CHAIN OF CUSTODY CHRONICLE

1. Shuttle Opened By: (print) Mark Darnas Date: 7/22/85 Time: 1100 AM
Signature: Mark Darnas Seal #: 37077 Intact:

I have received these materials in good condition from the above person.
2. Name: M. L. Danner Signature: M. L. Danner

Date: 7/22/85 Time: 1100 HR Remarks: Mark Danas

I have received these materials in good condition from the above person.
3. Name: _____ Signature: _____

Date: _____ Time: _____ Remarks: _____

4. Shuttle Sealed By: (print) Mark Danas Date: 7/23/85 Time: 10:30 AM
Signature: Mark Danas Seal #: 0031078 Intact:

ETC USE ONLY Opened By: Khanna Date: 7/24/85 Time: 4:00pm
Seal #: 37078 Condition: C+

FIELD PARAMETER FORM (CC2)

ETC JOB # J3284

Sample Point

W

B-19

Source Code

Sample Point ID

FIELD PROCEDURES

07/22/85

PURGE DATE
(YY MM DD)

09/30

START PURGE
(2400 hr Clock)

124

ELAPSED HRS

WATER VOL IN CASING
(Gallons)

1/10

VOLUME PURGED
(Gallons)

SAMPLING METHOD:

Bailer

Sampler Type



- A-Submersible Pump
B-ISCO
C-Bladder Pump

- D-Dipper/Bottle
E-Bailer
F-Scoop/Shovel

X-Other

(SPECIFY OTHER)

Sampler Material



- A-Teflon
B-Metal

- C-PVC
D-Plastic

X-Other

(SPECIFY OTHER)

Tubing Material



- A-Teflon
B-Tygon

- C-Polyethylene
D-Silicon

X-Other

(SPECIFY OTHER)

Sample Composited

Y/N

Procedure/Proportions

FIELD MEASUREMENTS

Well Elevation (ft/msl)

103501

Well Depth (ft)

13011

Depth to Ground water (ft)

120117

Sample Depth (non-well) (ft)

11111

Groundwater Elevation (ft msl)

1014814

1st pH (STD)1st spec. cond.um/cm³
at 25 °C (other parameter) value units2nd pH (STD)2nd spec. cond.um/cm³
at 25 °C (other parameter) value units3rd pH (STD)3rd spec. cond.um/cm³
at 25 °C (other parameter) value units4th pH (STD)4th spec. cond.um/cm³
at 25 °C (other parameter) value units Sample Temp (°C) Turbidity

NTU

FIELD COMMENTS

Sample Appearance: ~~soil~~ silty, brown color due to silt

Weather Conditions: Partly cloudy, 79°F

Other:

FILTERING: Use Chain of Custody (CC1) to indicate which bottles were filtered

Sampler: Mark Danas
(Print)Employer: CBS Records

I certify that sampling procedures were in accordance with applicable EPA state and corporate protocols.

Mark Danas 23 July 1985
(Signature)

**ETC ENVIRONMENTAL
TESTING AND CERTIFICATION**

CHAIN OF CUSTODY FORM (CC1)

Shuttle B

Seal No. 37075 ETC Job # 53284

Date Sealed 7/12/85 By: Vlance

Company: CBS Records

Attn.: mark donas

Facility/Site: _____

Phone: (404) 836-2343

Address: 5152 Columbia Drive

2020-2021 School Year | Page 10

Carrollton, GA 30117

SAMPLE IDENTIFICATION

Sample Point:	W-B-9	Your Sample Point ID	0723815	Start Date (YYMMDD)	1030	Start Time (2400 hr. clock)	0000	Elapsed Hours from sample.
Source Code (from below)	Nett weight	Sum Date (YYMMDD)						

Source Codes: Well ... (W) Outfall (O) Bottom Sediment (B) Surface Impoundment....(I) Leachate Collection Sys....(C) Other(X)
Soil ... (S) River/Stream ..(R) Generation Point(G) Treatment Facility(T) Lake/Ocean(L) Specify _____

SHUTTLE CONTENTS

CHAIN OF CUSTODY CHRONICLE

- | | | | | |
|--|--|-----------------------|----------------------------------|---------------------------|
| 1. | Shuttle Opened By: (print)
Signature: | Mark Danas | Date: 7/22/85
Seal #: 37075 | Time: 1100 HR
Intact: |
| I have received these materials in good condition from the above person. | | | | |
| 2. | Name:
Date: | Mark Danas
7/22/85 | Signature:
Remarks: | Mark Danas |
| I have received these materials in good condition from the above person. | | | | |
| 3. | Name:
Date: | | Signature:
Remarks: | |
| 4. | Shuttle Sealed By: (print)
Signature: | Mark Danas | Date: 7/23/85
Seal #: 0037076 | Time: 10:30 AM
Intact: |

ETC USE ONLY Opened By: Khanna Date: 7/24/85 Time: 11:00 am
Seal #: 37076 Condition: OK 4:00 pm

FIELD PARAMETER FORM (CC2)

FIELD PROCEDURES

07-22-85
PURGE DATE
(YY MM DD)0930
START PURGE
(2400 Hr. CLOCK)1241
ELAPSED HRSWATER VOL. IN CASING
Gallons1110
VOLUME PURGED
Gallons

SAMPLING METHOD: Bailer

Sampler Type	<input checked="" type="checkbox"/>	A-Submersible Pump B-ISCO C-Bladder Pump	D-Dipper/Bottle E-Bailer F-Scoop/Shovel	X-Other _____ <small>(SPECIFY OTHER)</small>
Sampler Material	<input checked="" type="checkbox"/>	A-Teflon B-Metal	C-PVC D-Plastic	X-Other _____ <small>(SPECIFY OTHER)</small>
Tubing Material	<input type="checkbox"/>	A-Teflon B-Tygon	C-Polyethylene D-Silicon	X-Other _____ <small>(SPECIFY OTHER)</small>
Sample Composited	<input type="checkbox"/> Y/N	Procedure/Proportions		

FIELD MEASUREMENTS

Well Elevation (ft/msl)

11013151011

Well Depth (ft)

13011

Depth to Ground water (ft)

20117

Sample Depth (non-well) (ft)

1111

Groundwater Elevation (ft msl)

10141819

1st	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	(STD)	1st	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	um/cm at 25°C			
2nd	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	(STD)	2nd	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	um/cm at 25°C			
3rd	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	(STD)	3rd	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	um/cm at 25°C			
4th	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	(STD)	4th	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	um/cm at 25°C			
	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	(°C)		<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	NTU			
	Sample Temp			Turbidity				
						(other parameter)	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	value
						(other parameter)	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	value
						(other parameter)	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	value
						(other parameter)	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	value
						(other parameter)	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	value
						(other parameter)	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	value

FIELD COMMENTS

Sample Appearance: Silty, brown color due to silt

Weather Conditions: Partly cloudy, 79°F

Other:

FILTERING: Use Chain of Custody (CC1) to indicate which bottles were filtered

Sampler: Mark Darnas
(Print)Employer: CBS Records

I certify that sampling procedures were in accordance with applicable EPA state and corporate protocols.

Mark Darnas 23 July 1985
(Date) (Signature)

LOG LINK 10496

LABORATORY CHRONICLE: Metals Department

Samples J 3284

	Chemist	Date
Hg Prep	<u>R Walls</u>	<u>8/6/85</u>
AA/ICAP Prep	<u>Joan Komarek</u>	<u>8/3/85</u>
Lab Supervisor	<u>Anthony W Macario</u>	Date <u>8-9-85</u>

JUL 31, 1985 11:54
SEQUENCE: 88195 ON CRN 15
CHANNEL 0

SUBSEQUENCE 1

METHOD
PCBB0

DIALG-PRG PARAM-FILE
/N

#WSHS #PMPS STOP
5, 5, 1

ISO POST-BTL# POST-#WSHS
NO , 0, 1

SAMPLES

	SAMPLE-NAME	BTL#	PROC-FILE	RAW-FILE	NDIL-F	STD-AMT	SMP-AMT
1	QC3331	,	1, PJ1595:FL	RJ1595:FL	100.00	1.0000	1.0000
2	AR1260-0.03	,	2, PJ1596:FL	RJ1596:FL	100000	1.0000	1.0000
3	AR1260-0.1	,	3, PJ1597:FL	RJ1597:FL	100000	1.0000	1.0000
4	QC33315	,	4, PJ1598:FL	RJ1598:FL	100.00	1.0000	1.0000
5	J0255	,	5, PJ1599:FL	RJ1599:FL	100.00	1.0000	1.0000
6	J02555	,	6, PJ1600:FL	RJ1600:FL	100.00	1.0000	1.0000
7	J0214	,	7, PJ1601:FL	RJ1601:FL	100.00	1.0000	1.0000
8	J0214R	,	8, PJ1602:FL	RJ1602:FL	100.00	1.0000	1.0000
9	J0304	,	9, PJ1603:FL	RJ1603:FL	200.00	1.0000	1.0000
10	J0305	,	10, PJ1604:FL	RJ1604:FL	200.00	1.0000	1.0000
11	J0306	,	11, PJ1605:FL	RJ1605:FL	200.00	1.0000	1.0000
12	J0307	,	12, PJ1606:FL	RJ1606:FL	200.00	1.0000	1.0000
13	J0312	,	13, PJ1607:FL	RJ1607:FL	200.00	1.0000	1.0000
14	AR1260-0.5	,	14, PJ1608:FL	RJ1608:FL	100000	1.0000	1.0000
15	AR1254-0.5	,	15, PJ1609:FL	RJ1609:FL	100000	1.0000	1.0000
16	AR1248-0.5	,	16, PJ1610:FL	RJ1610:FL	100000	1.0000	1.0000
17	AR1242-0.5	,	17, PJ1611:FL	RJ1611:FL	100000	1.0000	1.0000
18	J0191	,	18, PJ1612:FL	RJ1612:FL	100.00	1.0000	1.0000
19	J0197	,	19, PJ1613:FL	RJ1613:FL	100.00	1.0000	1.0000
20	J0193	,	20, PJ1614:FL	RJ1614:FL	189.00	1.0000	1.0000
21	J0316	,	21, PJ1615:FL	RJ1615:FL	200.00	1.0000	1.0000
22	J0258	,	22, PJ1616:FL	RJ1616:FL	100.00	1.0000	1.0000
23	J0274	,	23, PJ1617:FL	RJ1617:FL	100.00	1.0000	1.0000
24	AR1232-0.5	,	24, PJ1618:FL	RJ1618:FL	100000	1.0000	1.0000
25	AR1221-0.5	,	25, PJ1619:FL	RJ1619:FL	100000	1.0000	1.0000
26	AR1016-0.5	,	26, PJ1620:FL	RJ1620:FL	100000	1.0000	1.0000
27	AR1260-1.0	,	27, PJ1621:FL	RJ1621:FL	100000	1.0000	1.0000
28	<u>QC3335</u>	,	28, PJ1622:FL	RJ1622:FL	200.00	1.0000	1.0000
29	DMPST C	,	29, PJ1623:FL	RJ1623:FL	100000	1.0000	1.0000
30	AR1260-1.0	,	30, PJ1624:FL	RJ1624:FL	100000	1.0000	1.0000
31	QC33355	,	31, PJ1625:FL	RJ1625:FL	200.00	1.0000	1.0000
32	J1609	,	32, PJ1626:FL	RJ1626:FL	204.00	1.0000	1.0000
33	J16095	,	33, PJ1627:FL	RJ1627:FL	200.00	1.0000	1.0000
34	J1366	,	34, PJ1628:FL	RJ1628:FL	213.00	1.0000	1.0000
35	J1366R	,	35, PJ1629:FL	RJ1629:FL	200.00	1.0000	1.0000
36	J1615	,	36, PJ1630:FL	RJ1630:FL	200.00	1.0000	1.0000
37	J4012	,	37, PJ1631:FL	RJ1631:FL	200.00	1.0000	1.0000
38	J4017	,	38, PJ1632:FL	RJ1632:FL	200.00	1.0000	1.0000

J Stan 7/31/85

Kruehl 8/6/85

June G. Schaper 8/15/85

39	J1363	,	39,	PJ1633:FL	, RJ1633:FL	, 200.00,	1.0000,	1.0000
40	J1369	,	40,	PJ1634:FL	, RJ1634:FL	, 200.00,	1.0000,	1.0000
41	AR1260-0.5	,	41,	PJ1635:FL	, RJ1635:FL	, 100000,	1.0000,	1.0000
42	AR1254-0.5	,	42,	PJ1636:FL	, RJ1636:FL	, 100000,	1.0000,	1.0000
43	AR1248-0.5	,	43,	PJ1637:FL	, RJ1637:FL	, 100000,	1.0000,	1.0000
44	AR1242-0.5	,	44,	PJ1638:FL	, RJ1638:FL	, 100000,	1.0000,	1.0000
45	J1370	,	45,	PJ1639:FL	, RJ1639:FL	, 200.00,	1.0000,	1.0000
46	J1371	,	46,	PJ1640:FL	, RJ1640:FL	, 200.00,	1.0000,	1.0000
47	J1372	,	47,	PJ1641:FL	, RJ1641:FL	, 200.00,	1.0000,	1.0000
48	J2956	,	48,	PJ1642:FL	, RJ1642:FL	, 200.00,	1.0000,	1.0000
49	J3702	,	49,	PJ1643:FL	, RJ1643:FL	, 200.00,	1.0000,	1.0000
50	J3284	,	50,	PJ1644:FL	, RJ1644:FL	, 200.00,	1.0000,	1.0000
51	J2883	,	51,	PJ1645:FL	, RJ1645:FL	, 200.00,	1.0000,	1.0000
52	J4165	,	52,	PJ1646:FL	, RJ1646:FL	, 200.00,	1.0000,	1.0000
53	AR1232-0.5	,	53,	PJ1647:FL	, RJ1647:FL	, 100000,	1.0000,	1.0000
54	AR1221-0.5	,	54,	PJ1648:FL	, RJ1648:FL	, 100000,	1.0000,	1.0000
55	AR1016-0.5	,	55,	PJ1649:FL	, RJ1649:FL	, 100000,	1.0000,	1.0000
56	AR1260-5.0	,	56,	PJ1650:FL	, RJ1650:FL	, 100000,	1.0000,	1.0000
57	/E							

LABORATORY CHRONICLE: Sample Preparation Department

Sample Number	Log Link	Sample Vol. (ml)	Extract Vol. (ml)		Comments
			PEST	—	
J3285	10,409	500	5.0	/	AP 7
J3286	/	500	5.0	/	
J3287	/	500	5.0	/	
J3288	/	500	5.0	/	
J3289	/	500	5.0	/	
J3290	✓	500	5.0	/	↓
J2956	10,441	500	5.0	/	AP8
J3702	10,470	460	5.0	/	
J3284	10,496	480	5.0	/	
J2883	10,484	500	5.0	/	
J3096	10,524	500	5.0	/	
J3097	/	500	5.0	/	
J3098	/	500	5.0	/	
J3099	/	500	5.0	/	
J3102	/	500	5.0	/	
J3101	/	500	5.0	/	
J3103	/	500	5.0	/	
J3104	✓	460	5.0	/	↓
QC 3413	(S3)	500	5.0	/	
QC 3413	(S2)	500	5.0	/	
J3101	(S2)	450	5.0	/	
QC 3413		500	5.0	/	
QC 3413	(S1)	500	5.0	/	
J3102	(S2)	450	5.0	/	
J3103	(S2)	450	5.0	/	
J3288	R	500	5.0	/	

FRACTION	SPIKE		
	Amt (ml)	Conc.	Lot #
(S1) ARS-Pest (EC) SPIKE	0.1	14025 µg/ml	11,044
		18250 µg/ml	
(S2) DR8-Pest (EC) FICHEMICAL	0.1	2025 µg/ml	10,636
		10100 µg/ml	
(S3) AP8-Pest (FPD) SPIKE	0.2	4012.5 µg/ml	10,640
		20625 µg/ml	
		10125 µg/ml	

QC Batch # 3413

Analysis AP7/PEST(EC)

AP8/PEST(EC)

AP8/PEST(FPD)

Matrix H₂O

Turnaround Normal

Date 7/30/85

Extraction Method:

Sep. Funnel ✓

Continuous

Soxhlet

Other

COMMENTS

Set-up: 9 Aug 85 / JF 7/30/85	UPD/Supervisor: Dan L. 7/31/85-
Spec.: 1/1 mm 7/31/85	Spike/Surr. Verified: Lee Burkhardt 7/31/85

AUG 5, 1985 14:28
SEQUENCE: FP3413 ON CRN 15
CHANNEL 8 QC 3413
SUBSEQUENCE 1 Pts: FPD
 primary

a. Miggiaiu 8/5/85
FPD: a. Miggiaiu 8/6/85

METHOD
AP8FPD

DIALG-PRG PARAM-FILE
/N

#WSHS #PMPS STOP
5, 5, 1

ISO POST-BTL# POST-#WSHS
NO , 0, 1

SAMPLES

	SAMPLE-NAME	BTL#	PROC-FILE	RAW-FILE	ZDIL-F	STD-AMT	SMP-AMT
1	HEXANE	,	1, PA2655:FW,	RA2655:FW,	100.00,	1.0000,	1.0000
2	QC3413	,	2, PA2656:FW,	RA2656:FW,	1000.0,	1.0000,	1.0000
3	AP8FPD-B	,	3, PA2657:FW,	RA2657:FW,	100000,	1.0000,	1.0000
4	QC3413S3	,	4, PA2658:FW,	RA2658:FW,	1000.0,	1.0000,	1.0000
5	J3101	,	5, J3101S:FW,	J3101T:FW,	1000.0,	1.0000,	1.0000
6	J3101S3	,	6, PA2659:FW,	RA2659:FW,	1111.0,	1.0000,	1.0000
7	J3288	,	7, J3288S:FW,	J3288T:FW,	1000.0,	1.0000,	1.0000
8	J3288R	,	8, PA2660:FW,	RA2660:FW,	1000.0,	1.0000,	1.0000
9	J2956	,	9, J2956S:FW,	J2956T:FW,	1000.0,	1.0000,	1.0000
10	J3702	,	10, J3702S:FW,	J3702T:FW,	1087.0,	1.0000,	1.0000
11	J3284	,	11, J3284S:FW,	J3284T:FW,	1042.0,	1.0000,	1.0000
12	J2883	,	12, J2883S:FW,	J2883T:FW,	1000.0,	1.0000,	1.0000
13	J3096	,	13, J3096S:FW,	J3096T:FW,	1000.0,	1.0000,	1.0000
14	AP8FPD-A	,	14, PA2661:FW,	RA2661:FW,	100000,	1.0000,	1.0000
15	J3097	,	15, J3097S:FW,	J3097T:FW,	1000.0,	1.0000,	1.0000
16	J3098	,	16, J3098S:FW,	J3098T:FW,	1000.0,	1.0000,	1.0000
17	J3099	,	17, J3099S:FW,	J3099T:FW,	1000.0,	1.0000,	1.0000
18	J3102	,	18, J3102S:FW,	J3102T:FW,	1000.0,	1.0000,	1.0000
19	J3103	,	19, J3103S:FW,	J3103T:FW,	1000.0,	1.0000,	1.0000
20	J3104	,	20, J3104S:FW,	J3104T:FW,	1087.0,	1.0000,	1.0000
21	AP8FPD-C	,	21, PA2662:FW,	RA2662:FW,	100000,	1.0000,	1.0000
22	/E						

R. Baker
8-11-85

AUG 6, 1985 13:58
SEQUENCE: S3413P ON CRN 15
CHANNEL 8

a.m. 8/6/85

Pest(EC) A. megaceia 8/6/85

SUBSEQUENCE 1

METHOD
AP8PST

DIALG-PRG PARAM-FILE
/N

#WSHS #PMPS STOP
5, 5, 1

ISO POST-BTL# POST-#WSHS
NO, 0, 1

R. Baker
8-14-85

SAMPLES

	SAMPLE-NAME	BTL#	PROC-FILE	RAW-FILE	XDIL-F	STD-AMT	SMP-AMT
1	HEXANE	,	1, PA2670:GB,	RA2670:GB,	100.00,	1.0000,	1.0000
2	QC3413	,	2, PA2671:GB,	RA2671:GB,	1000.0,	1.0000,	1.0000
3	AP8PST-B	,	3, PA2672:GB,	RA2672:GB,	100000,	1.0000,	1.0000
4	QC3413S1	,	4, PA2673:GB,	RA2673:GB,	1000.0,	1.0000,	1.0000
5	J3103	,	5, PA2674:GB,	RA2674:GB,	1000.0,	1.0000,	1.0000
6	J3103S1	,	6, PA2675:GB,	RA2675:GB,	1111.0,	1.0000,	1.0000
7	J3288	,	7, PA2676:GB,	RA2676:GB,	1000.0,	1.0000,	1.0000
8	J3288R	,	8, PA2677:GB,	RA2677:GB,	1000.0,	1.0000,	1.0000
9	J2956	,	9, PA2678:GB,	RA2678:GB,	1000.0,	1.0000,	1.0000
10	J3702	,	10, PA2679:GB,	RA2679:GB,	1087.0,	1.0000,	1.0000
11	J3284	,	11, PA2680:GB,	RA2680:GB,	1042.0,	1.0000,	1.0000
12	J2883	,	12, PA2681:GB,	RA2681:GB,	1000.0,	1.0000,	1.0000
13	J3096	,	13, PA2682:GB,	RA2682:GB,	1000.0,	1.0000,	1.0000
14	AP8PST-A	,	14, PA2683:GB,	RA2683:GB,	100000,	1.0000,	1.0000
15	J3097	,	15, PA2684:GB,	RA2684:GB,	1000.0,	1.0000,	1.0000
16	J3098	,	16, PA2685:GB,	RA2685:GB,	1000.0,	1.0000,	1.0000
17	J3099	,	17, PA2686:GB,	RA2686:GB,	1000.0,	1.0000,	1.0000
18	J3102	,	18, PA2687:GB,	RA2687:GB,	1000.0,	1.0000,	1.0000
19	J3101	,	19, PA2688:GB,	RA2688:GB,	1000.0,	1.0000,	1.0000
20	J3104	,	20, PA2689:GB,	RA2689:GB,	1087.0,	1.0000,	1.0000
21	J3285	,	21, PA2690:GB,	RA2690:GB,	1000.0,	1.0000,	1.0000
22	J3286	,	22, PA2691:GB,	RA2691:GB,	1000.0,	1.0000,	1.0000
23	J3287	,	23, PA2692:GB,	RA2692:GB,	1000.0,	1.0000,	1.0000
24	J3289	,	24, PA2693:GB,	RA2693:GB,	1000.0,	1.0000,	1.0000
25	J3290	,	25, PA2694:GB,	RA2694:GB,	1000.0,	1.0000,	1.0000
26	AP8PST-C	,	26, PA2695:GB,	RA2695:GB,	100000,	1.0000,	1.0000
27	HEXANE	,	27, PA2696:GB,	RA2696:GB,	100.00,	1.0000,	1.0000
28	AP8ISO-B	,	28, PA2697:GB,	RA2697:GB,	100000,	1.0000,	1.0000
29	QC3413S2	,	29, PA2698:GB,	RA2698:GB,	1000.0,	1.0000,	1.0000
30	HEXANE	,	30, PA2699:GB,	RA2699:GB,	100.00,	1.0000,	1.0000
31	J3102S2	,	31, PA2700:GB,	RA2700:GB,	1111.0,	1.0000,	1.0000
32	HEXANE	,	32, PA2701:GB,	RA2701:GB,	100.00,	1.0000,	1.0000
33	AP8ISO-A	,	33, PA2702:GB,	RA2702:GB,	100000,	1.0000,	1.0000
34	AP8ISO-C	,	34, PA2703:GB,	RA2703:GB,	100000,	1.0000,	1.0000
35	/E						

AUG 7, 1985 14:58
SEQUENCE: D3413P ON CRN 16
CHANNEL 8

SUBSEQUENCE 1

Pest(EC)
Dilutions

A. Mufacciu 8/7/85+
8/8/85

METHOD
AP8PST

DIALG-PRG PARAM-FILE
/N

#WSHS #PMPS STOP
5, 5, 1

ISO POST-BTL# POST-#WSHS
NO , 0, 1

J. Baker
8/14/85

SAMPLES

	SAMPLE-NAME	BTL#	PROC-FILE	RAW-FILE	XDIL-F	STD-AMT	SMP-AMT
1	AP8PST-B	,	1, PA2730:GJ,	RA2730:GJ,	100000,	1.0000,	1.0000
2	J2883 1:100	,	2, PA2731:GJ,	RA2731:GJ,	100000,	1.0000,	1.0000
3	J3096 1:25	,	3, PA2732:GJ,	RA2732:GJ,	25000.,	1.0000,	1.0000
4	J3097 1:100	,	4, PA2733:GJ,	RA2733:GJ,	100000,	1.0000,	1.0000
5	J3098 1:100	,	5, PA2734:GJ,	RA2734:GJ,	100000,	1.0000,	1.0000
6	J3099 1:50	,	6, PA2735:GJ,	RA2735:GJ,	50000.,	1.0000,	1.0000
7	AP8PST-A	,	7, PA2736:GJ,	RA2736:GJ,	100000,	1.0000,	1.0000
8	/E						

LABORATORY CHRONICLE: GC-MS Department

DATE 850801 SHIFT _____
FRACTION VOA
INSTRUMENT D
TUNE FILE ~~WOBADS~~ APE 102
SEQUENCE FILE _____
METHOD FILE VOAD
ID FILE DAPR8PT
ANALYST(S) W. Debaa
SUPERVISOR _____
BATCH #'s Q13555

(PLEASE INITIAL)

CURRENT CS05 STATUS	STANDARDS UPDATED
ACQ	DATE
DC	850802 JL
WIP	BY

NAME	DATA FILE	uL INJ	ALS #	DIL	TAPE #	SPECIALS (WRITE A-TYPE)	PLUS Y/N
P-BFB	>D0071	1ml				OK	
QC 3555V	>D0072	5ml	1				
QC 3555VS	>D0073		2			5ml ABC/Sewer + AP8	
QC 3555VS	>D0074		3			10ml ABC/Sewer + AP8	
QC 3555VS	>D0075	↓	4			30ml ABC/Sewer + AP8	
P-BFB	>D0076					119 - 140	
J3098V5	>D0077		1				
J3098V	>D0078		2				
J3702V	>D0079		3				
J3702VR	>D0080		4				
J3284V	>D0081		5				
J2883V	>D0082		6	1:10		Rerun 1:100	
J3100V	>D0083	↓	7	1:10		Rerun for carry over	
J4060V	>D0084	↓	8	1:1000		Rerun 1:50	
(P) J3800V	>D0085			1:5		Rerun straight	
(P) J2883V	>D0086			1:100		OK	
P-BFB	>D0087					OK (122-112) (0853 mts)	
QC 3555VS	>D0088	5ml	1			OK 5ml ABC/Sewer + AP8.	
(P) J4060V	>D0089	100ml	1:50 ^{PS}	1:50		OK	
(P) J3100V	>D0090	5ml	2			OK	
J2258V	>D0091		3			OK	
J2257V	>D0092		4			OK	
J3097V	>D0093	↓	5			Rerun 1:10 for tailout.	
(P) J3097V	>D0094	.5ml	1	1:10		OK	

Sample Number	Log Link	Sample Vol. (ml)	Extract	Comments
			Vol. (ml)	
J2257	10,316	500	2.5	
J2258	↓		2.5	
J2956	10,441		2.5	
J3702	10,470		2.5	
J3284	10,496		2.5	
J2883	10,484		2.5	
J3046	10,524		2.5	
J3097			5.0	SAMPLE VISCOUS @ 2.5 ml F.V.
J3098			2.5	
J3099			2.5	
J3101			2.5	
J3102			2.5	
J3103			2.5	
J3104	↓		2.5	
QC 3401		500	2.5	
QC 3401	S		2.5	
J3104	S		2.5	
J3102	R	↓	2.5	

FRACTION	SPIKE		
	Amt (ml)	Conc.	Lot #
2,4-D	0.5	100µg/ml	11)114
2,4,5-TP		10µg/ml	
2,4,5-T	0.5	10µg/ml	11)115

QC Batch #	3401
Analysis	AP8/HGRB (EC)
Matrix	H ₂ O
Turnaround	Normal
Date	7/29/85
Extraction Method:	8/26
Sep. Funnel	
Continuous	
Soxhlet	
Other	
COMMENTS	

Blend

	SURROGATE		
	Amt. (ml)	Conc.	Lot #

Set-up: Wendy Diangeloni 7/29/85
 Conc.: Wendy Diangeloni 7/29/85
 ETC FORM #105

UPD/Supervisor: Diane - 9/26/85
 Spike/Surr. Verified: Karen Albrecht 7/29/85

AUG 13, 1985 14:22

SEQUENCE: AP8HER ON CRN 16

CHANNEL #

SUBSEQUENCE 1

QC 3401

3458
3423

A. mucifacium 8/13/85

} AP8Her primary

METHOD

HERB

DIALG-PRG PARAM-FILE

/N

#WSHS #PMPS STOP

5, 5, 1

ISO POST-BTL# POST-#WSHS

NO , 0, 1

SAMPLES

	SAMPLE-NAME	BTL#	PROC-FILE	RAW-FILE	NDIL-F	STD-AMT	SMP-AMT
1	HEXANE	,	1, PA2825:HA,	RA2825:HA,	100.00,	1.0000,	1.0000
2	QC3401	,	2, PA2826:HA,	RA2826:HA,	500.00,	1.0000,	1.0000
3	HERB- EAB	,	3, PA2827:HA,	RA2827:HA,	100000,	1.0000,	1.0000
4	245T- 5.0 10	,	4, PA2828:HA,	RA2828:HA,	100000,	1.0000,	1.0000
5	QC3401S	,	5, PA2829:HA,	RA2829:HA,	500.00,	1.0000,	1.0000
6	J3104	,	6, PA2830:HA,	RA2830:HA,	500.00,	1.0000,	1.0000
7	J3104S	,	7, PA2831:HA,	RA2831:HA,	500.00,	1.0000,	1.0000
8	J3102	,	8, PA2832:HA,	RA2832:HA,	500.00,	1.0000,	1.0000
9	J3102R	,	9, PA2833:HA,	RA2833:HA,	500.00,	1.0000,	1.0000
10	J2257	,	10, PA2834:HA,	RA2834:HA,	500.00,	1.0000,	1.0000
11	J2258	,	11, PA2835:HA,	RA2835:HA,	500.00,	1.0000,	1.0000
12	J2956	,	12, PA2836:HA,	RA2836:HA,	500.00,	1.0000,	1.0000
13	J3702	,	13, PA2837:HA,	RA2837:HA,	500.00,	1.0000,	1.0000
14	J3224	,	14, PA2838:HA,	RA2838:HA,	500.00,	1.0000,	1.0000
15	HERB- PBA	,	15, PA2839:HA,	RA2839:HA,	100000,	1.0000,	1.0000
16	245T- 5.0 10	,	16, PA2840:HA,	RA2840:HA,	100000,	1.0000,	1.0000
17	J2883	,	17, PA2841:HA,	RA2841:HA,	500.00,	1.0000,	1.0000
18	J3096	,	18, PA2842:HA,	RA2842:HA,	500.00,	1.0000,	1.0000
19	J3097	,	19, PA2843:HA,	RA2843:HA,	500.00,	1.0000,	1.0000
20	J3098	,	20, PA2844:HA,	RA2844:HA,	500.00,	1.0000,	1.0000
21	J3099	,	21, PA2845:HA,	RA2845:HA,	500.00,	1.0000,	1.0000
22	J3101	,	22, PA2846:HA,	RA2846:HA,	500.00,	1.0000,	1.0000
23	J3103	,	23, PA2847:HA,	RA2847:HA,	500.00,	1.0000,	1.0000
24	HERB-C	,	24, PA2848:HA,	RA2848:HA,	100000,	1.0000,	1.0000
25	245T-10,0	,	25, PA2849:HA,	RA2849:HA,	100000,	1.0000,	1.0000

K. Banks
K. 8/19/85

Set-up: Wanda Elmore 8/11/85

Conc.: Wards Elementary 8/11/85

UPD/Supervisor: Dan Feltman 8/12/25

Spike/Surr. Verified:

QC Batch #	3Y29
Analysis	<u>GAP8/FORM</u>
Matrix	<u>H₂O</u>
Turnaround	<u>Normal</u>
Date	<u>8/10</u>
<u>Extraction Method:</u>	
Sep. Funnel	<u> </u>
Continuous	<u> </u>
Soxhlet	<u> </u>
Other	<u>✓</u>
COMMENTS	

SURROGATE	
Conc.	Lot #

LABORATORY CHRONICLE: GC-MS Department

DATE August 17th 85 SHIFT _____
FRACTION Douvr. / Form.
INSTRUMENT I
TUNE FILE MTI 001
SEQUENCE FILE KV15.B
METHOD FILE BNPI
ID FILE FORM
ANALYST(S) A. Isakson

SUPERVISOR P.T. - b
BATCH #'s QB 3429 - Derry

(PLEASE INITIAL)

CURRENT CS05 STATUS	STANDARDS UPDATED
ACQ	DATE
WIP	BY

CB FORM 14 CB-444

FRACTION	SPIKE		
	Amt (ml)	Conc.	Lot #

FLUORANTHENE	1.0	0.084 μM	10.857
PARATHION	1.0	0.4 μg/ml	10.858
Z-45	1.0	0.4 μg/ml	10.859

(vishnugopha)

Set-up: Ken Bapkin 8/9/85

Conc.: Ken Brink 8/10/85

UPD/Supervisor: Karen Alvesten 8/10/05

Spike/Surr. Verified: Bombo & 35. Elians 8/9/85

QC Batch # 3424

Analysis CSAP8 / 8/6/10

Matrix H₂O

Turnaround NORMAL

Date 8/8/85

Extraction Method:

Sep. Funnel _____

Continuous _____

Soxhlet _____

Other sep - sep R

COMMENTS

Ken did not
designate on
labels which
J 9284 was the
SPIKE. (MA) 8/8/85

Also, he did not
differentiate
between the
QC + QCS. Was
repeated 8/10/85.

Analyst: Paula E. Koleno
8/23/85.

SURROGATE		
ml)	Conc.	Lot #

LABORATORY CHRONICLE: Sample Preparation Department

QC Batch # 3427

Analysis GAP8 18630

Matrix H₂O

Turnaround NORMAL

Date 2/13/85

Extraction Method:

Sep. Funnel _____

Continuous _____

Soxhlet _____

Other _____ ✓

COMMENTS

FRACTION	SPIKE		
	Amt (ml)	Conc.	Lot #
PHENOL	100 uL	40 v/v/dl	10,594

				SURROGATE		
				Amt. (ml)	Conc.	Lot #

Set-up: K. Borgia 8/13/85

UPD/Supervisor: Don Warr 8/16/85

Conc.: K, ~~Fe~~, 8/19/85

Spike/Surr. Verified: Dan Wootton 8/18/05

LABORATORY CHRONICLE: Sample Preparation Department

Set-up: K. Bachman 8/20/85
- L. Bachman 8/21/85

UPD/Supervisor:

Editor/Cust. Verified: D. J. L. 8/22/95

Paul R. Horner 8/21/85

~~Q-1~~ 1/2c/65

LABORATORY CHRONICLE: GC-MS Department

DATE 950808/9 SHIFT 1
 FRACTION A9748 BNA
 INSTRUMENT E
 TUNE FILE TME MT EOD1
 SEQUENCE FILE TME
 METHOD FILE BNDE
 ID FILE EAP8E, FBNA8
 ANALYST(S) Thomas Mancuso
 SUPERVISOR Pearl Frank
 BATCH #'S QC3418, QC3450

(PLEASE INITIAL)

CURRENT CS05 STATUS		STANDARDS UPDATED	
ACQ	<u>TM</u>	DATE	<u>850810</u>
WIP		BY	<u>Thomas Mancuso</u>

CBAPIB, CBAPIA, CIBIUM8

STANDARD	CONC PPM	LOT NO.	LOT VOL.
ISTD	400	11131	100/5ml
DFTPP	25	10941	2ml
BNCalib III	200	11026	1ml
" II	100	11151	/
" I	60	11239	/
AP8BN1 High	300	9899	/
" Med	100	9894	/
" Low	30	9895	/
AP8BN1 High	500	9891	/
" Med	300	9896	/
" Low	100	9897	/
Calib AP8 Acid	300	10856	/
" Med	150/100	10855	/
" Low	60	10854	/
AP8 Acid Calib	300	9893	/
" II	100	9899	/
" I	30	9899	/

NAME	DATA FILE #	uL INJ	ALS #	DIL	TAPE #	SPECIALS (WRITE A-TYPE)	PLUS Y/N
DFTPP	DE0696		1				
BNCalib III	0697		2				
" " II	0698		3				
" " I	0699		4				
AP8BN1 High	0700		5				
" Med	0701		6				
" Low	0702		7				
AP8BN1 High	0703		8				
" Med	0704		9				
" Low	0705		10				
Acid Calib III	0706		11				
" " II	0707		12				
" " I	0708		13				
AP8 Acid III	0709		14				
" II	0710		15				
" I	0711		16				
DFTPP	0712		17				
QC3418C	0713		18	-		Acid	
QC3418CS	0714		19	-			
J3289CS	0715		20	-			
J2257C	0716		21	-			
J2258C	0717		22	-			
J3285C	0718		23	.			
J3286C	0719		24	-			
J3287C	0720		25	-			

LABORATORY CHRONICLE: GC-MS Department

DATE 950809 SHIFT _____
FRACTION _____
INSTRUMENT _____
TUNE FILE _____
SEQUENCE FILE _____
METHOD FILE _____
ID FILE _____
ANALYST(S) Thomas/Moncur
SUPERVISOR Dee Trant
BATCH #'s _____

CURRENT CS05 STATUS	STANDARDS UPDATED	
ACQ	24	DATE
WIP		BY

NAME	DATA FILE E	UL INJ	ALS #	DIL	TAPE #	SPECIALS (WRITE A-TYPE)	PLUS Y/N
J3288C	>E0721		26	5		Acid	
J3288CR	0722		27				
J3289C	0723		28	-			
J3290C	0724		29	1			
J2956C	0725		30	1			
J3702C	0726		31	1			
J3284C	0727		32	1			
OF TPP	0728		33	1			
Acid Calb II	0729		34	1			
BN Calb II	0730		35	1			
QB3418C	0731		36	-		BN	
QB3418CS	0732		37	1			
J3289CS	0733		38	-			
J2257C	0734		39	1			
J2258C	0735		40	1			
J3285C	0736		41	1			
J3286C	0737		42	1			
J3287C	0738		43	1			
J3288C	0739		44	1			
J3289C	0740		45	1			
J3290C	0741		46	1			
J3288CB	0742		47	1			
J3702C	0743		48	1			
J2956C	0744		49	1			
J3702C	0745		50	1			

LABORATORY CHRONICLE: GC-MS Department

DATE 850811/12 SHIFT _____
FRACTION AP8 DER
INSTRUMENT J
TUNE FILE MTJ001
SEQUENCE FILE TMJ
METHOD FILE BNPJ1
ID FILE AP8 DER
ANALYST(S) ZM/jmcnew

SUPERVISOR
BATCH #'s Q03418, Q03430, Q03353

(PLEASE INITIAL)

CURRENT CS05 STATUS		STANDARDS UPDATED	
ACQ	DA	DATE	
WIP		BY	

NAME	DATA FILE	uL INJ	ALS #	DIL	TAPE #	SPECIALS (WRITE A-TYPE)	PLUS Y/N
DFTPP	754687		1				
Hexachlorophene deriv 300 ppm	754688		2				
" 150 ppm	4689		3				
" 50 ppm	4690		4				
QC3418C	4691		5				
QC3418CS	4692		6				
J3289CS	4693		7				
J2257C ⁽⁵⁰⁰⁾	4694		8				
J2258C ⁽⁵⁰⁰⁾	4695		9				
J3285C ⁽⁵⁰⁰⁾	4696		10				
J3286C	4697		11				
J3287C	4698		12				
J3288C	4699		13				
J3288CR	4700		14				
J3289C	4701		15				
J3290C	4702		16				
J2956C	4703		17				
DFTPP	4704		18				
Hexachlorophene 150 ppm STO	4705		19				
J3702C	4706		20				
J3284C	4707		21				
QC3353C	4708		22				
J4315CS	4709		23				
J4315C	4710		24				
J3701	4711		25				

TITLE: Cyanides

BOOK NO. CV 7

PAGE NO. 29

SAMPLE #	CONE.
1 J29510	ND
2 J3102	ND
3 J3284	ND
4 J2883	BMDL
5 J3096	BMDL
6 J3097	.29 ^{1:2}
7 J3098	BMDL
8 J3099	BMDL
9 J3100	.0491
10 J4040	7.13 ^{1:100}
11 J3640	2.96 ^{1:10}
12 J3647	81.0 ^{1:500}
13 J3653	97.9 ^{1:1000}
14 J3655	5.45 ^{1:50}
15 J3656	.287 ^{1:100}
16 J3658	9.47 ^{1:100}
17 J4759	46.3 ^{1:100}
18 J3644	35.8 ^{1:100}
19 J3646	1.84 ^{1:10}
20 J3659	1.73 ^{1:10}
21 J4059	9.15 ^{1:50}
22	
23	
24	
25	
26	
27	
28	
29	

SAMPLE #	CONE.
30	
31	
32	
33	
34	
35	
36	
37	
38	
39	
40	
41	
42	
43	
44	
45	
46	
47	
48	
49	
50	
51	
52	
53	
54	
55	
56	
57	
58	

UNITS: MG/L

Detection limit = 0.025 mg/l

COMMENTS:

< = less than 0.025 mg/l

Blank Recovery = .195/.200 = 98%

EPA Std. Rec. = .197/.224 = 88%

SPIKE RECOVERY

SAMPLE #: J3100 1:200

REPLICATE 1: .0473

REPLICATE 2: .0510

MEAN : .0491

STD. DEV. :

SPIKE VALUE: 0.200

REPLICATE 1: 0.199

REPLICATE 2:

MEAN :

STD. DEV. :

S. RECOVERY :

99%

REPLICATE DATA:

SAMPLE #: _____

REPLICATE 1: _____

REPLICATE 2: _____

MEAN :

STD. DEV. :

ANALYST: Deborah Kay
DATE : 8/29/83

VERIFIED BY: [Signature]

LABORATORY CHRONICLE: Sample Preparation Department

Sample Number	Log Link	Sample Vol. (ml)	Extract Vol. (ml)	Comments
J2257	10,316			ALIQUOTED
J2258	↓		1	
J3285	10,409		1	
J3286	↓	1		
J3287	↓	1		
J3288	↓	1		
J3289	↓	1		
J3290	↓	1		
J2956	10,441	✓		
J3702	10,470	✓		
J2883	10,484	✓		
J3284	10,496	✓		
J3096	10,524	✓		
J3097	↓	1		
J3098	↓	1		
J3099	↓	1		
J3101	↓	1		
J3102	↓	1		
J3103	↓	1		
J3104	↓	1		
J3097	(15)	1	↓	
QC 3382	1			ALIQUOTED
QC 3382	1 S			
J3098	1 S		↓	
J3096	1 R		↓	

FRACTION	SPIKE		
	Amt (ml)	Conc.	Lot #
AP8-HPLC DAT	1.0 ml	8025ug/ml	11,020
Spike I			
AP8-HPLC DAT	1.0 ml	2025ug/ml	11,015
Spike II			

QC Batch # 3282

Analysis AP8/DAILC 210
AP8/DAILC 254
AP7/DAILC

Matrix H₂O

Turnaround Normal

Date 7/29/85

Extraction Method:

Sep. Funnel _____

Continuous _____

Soxhlet _____

Other Aliquot

COMMENTS

* Use BRAVO New Spikes Only.

- Aliquot ≈ 2ml into vial — label vial with "Sharpie" — not with tape labels.

- QC is hydropure H₂O only.

- MS & QCS: Add 1.0 ml of both AP8/DAILC stds to 5.0 ml volumetric. Bring to 5.0 ml with sample & hydropure H₂O respectively.

- Mix; aliquot into vial.

Analyst Name: *[Signature]*
Date: 9/12 9/20-23/85

SURROGATE		
Amt. (ml)	Conc.	Lot #

Set-up: 5 Meister 7/29/85

UPD/Supervisor: 5 Meister 7/29/85

LABORATORY CHRONICLE: Sample Preparation Department

Sample Number	Log Link	Sample Vol. (ml)	Extract Vol. (ml)	Comments
J3702	10,470	1000	1.0	
J3284	10,496	1000	1.0	
J2883	10,484	1000	1.0	
J3096	10,524	440	1.0	
J3097		930	1.0	
J3098		940	1.0	
J3099		830	1.0	
J3101		900	1.0	
J3102		940	1.0	
J3103		950	1.0	
J3104	✓	870	1.0	
J3100	10,557	440	1.0	
X J4995	10,625	run as soil MB 8/7/85		
C J3649	10,668	890	10ml	
C J3651		970	10ml	
C J3657	✓	800	1.0	
QC 3457		1000	1.0	
QC 3457 S		1000	1.0	
J3104 S		950	1.0	
J3103 R		980	1.0	

FRACTION	SPIKE		
	Amt (ml)	Conc.	Lot #
APR EXT LC I	0.5	1000 µg/ml	10,897
APR EXT LC II	0.5	1000 µg/ml	10,898
APR EXT LC III	0.5	1000 µg/ml	10,899
APR EXT TLC IV	0.5	1000 µg/ml	10,900

QC Batch # 3457

Analysis AP8/EXT LC

Matrix H₂O

Turnaround Normal

Date 3/8/85 - 8/7/85 MB

Extraction Method:

Sep. Funnel ✓

Continuous ✓

Soxhlet

Other

COMMENTS

X Use the sample designated for turnaround.
There is 100ml of sample w/ 9.00 ml H₂O

Run as soil MB 8/7/85

analyzed 9/26-27/85 /HLL

SURROGATE		
Amt. (ml)	Conc.	Lot #

Set-up: John Tabbate 8/8/85 Initial 8/8/85 UPD/Supervisor: Diana Gaudin 8/3/85

Conc.: 0.016 C. 4/1/85

LABORATORY CHRONICLE: GC-MS Department

DATE 850926 SHIFT _____
FRACTION V09
INSTRUMENT D
TUNE FILE APB10Z
SEQUENCE FILE _____
METHOD FILE APB HPT
ID FILE D APB DI
ANALYST(S) Ralph Quinn

SUPERVISOR [Signature]
BATCH #'s QV 3774

(PLEASE INITIAL)

CURRENT CS05 STATUS	STANDARDS UPDATED
ACQ	DATE
WIP	BY