

5-4-92

-10

DOC930027

FINAL

DATA EVALUATION REPORT

OIL OF CITRONELLA

Study Type: Mutagenicity: Salmonella typhimurium/Mammalian Microsome  
Mutagenicity Assay

Prepared for:

Health Effects Division  
Office of Pesticide Programs  
Environmental Protection Agency  
1921 Jefferson Davis Highway  
Arlington, VA 22202

Prepared by:

Clement International Corporation  
9300 Lee Highway  
Fairfax, VA 22031-1207

Principal Reviewer Lynne T. Haber Date 5/4/92  
Lynne T. Haber, Ph.D.  
Independent Reviewer Sharon A. Segal for Date 5/4/92  
Nancy E. McCarroll, B.S.  
QA/QC Manager Sharon A. Segal Date 5/4/92  
Sharon Segal, Ph.D.

Contract Number: 68D10075  
Work Assignment Number: 1-76  
Clement Number: 91-247  
Project Officer: James Scott

GUIDELINE SERIES 84: MUTAGENICITY  
SALMONELLA

EPA Reviewer: J. Thomas McClintock, Ph.D.  
Reregistration Section, Science Analysis and  
Coordination Branch (H7509C)  
EPA Section Head: Albin Kocialski, Ph.D.  
Reregistration Section, Science Analysis and  
Coordination Branch (H7509C)

Signature: J. Thomas McClintock  
Date: 5/5/92  
Signature: A. Kocialski  
Date: 5/8/92

DATA EVALUATION REPORT

STUDY TYPE: Mutagenicity: Salmonella typhimurium/mammalian microsome  
mutagenicity assay

EPA IDENTIFICATION Numbers:

Tox Chem. Number: 21901

MRID Number: 421513-10

TEST MATERIAL: TREO SPF 15 lotion

SYNONYMS: Oil of citronella

SPONSOR: Primavera Laboratories, Inc., New York, NY

STUDY NUMBER: 063629-2

TESTING FACILITY: United States Testing Company, Inc., Hoboken, NJ

TITLE OF REPORT: Ames Salmonella Microsome Mutagenesis Assay on TREO SPF 15  
Lotion

AUTHOR: Wang, X.M.

REPORT ISSUED: December 2, 1991

CONCLUSIONS--EXECUTIVE SUMMARY: No conclusions can be reached from the  
Salmonella typhimurium/mammalian microsome plate incorporation assay conducted  
with TREO SPF 15 lotion. The study was seriously compromised for the  
following reasons:

- 1--Unacceptably high spontaneous reversion counts for strains TA98 and  
TA100.
- 2--Marked reductions in the sensitivity of strains TA98 and TA100 to  
direct-acting mutagens.
- 3--The use of an excessive concentration of the S9-activated positive  
control (10 µg/plate 2-anthramine).

# SALMONELLA

(See Section C, Reported Results, for a detailed discussion.)

Based on the above considerations, it was concluded that the study was unacceptable, and, therefore, did not satisfy Guideline requirements (§84-2) for genetic effects Category I, Gene Mutations.

STUDY CLASSIFICATION: Unacceptable. The study should be repeated using established procedures for the Salmonella typhimurium/mammalian microsome plate incorporation assay.<sup>1</sup> The purity of the test compound should also be included with the report of the repeat study.

## A. MATERIALS:

### 1. Test Material: TREO SPF 15 lotion

Description: Off-white emulsion  
 Lot number: Not reported  
 Purity: Not reported  
 Receipt date: Not reported  
 Stability: Not reported  
 Contaminants: None listed  
 Solvent used: Dimethyl sulfoxide (DMSO)  
 Other provided information: The storage conditions and frequency of dose solution preparation were not reported.

### 2. Control Materials:

Negative: Not done

Solvent/final concentration: DMSO/≤0.1 mL

Positive:

Nonactivation:

N-Methyl-N-nitro-	<u>10</u>	µg/plate TA1535, TA100
N-nitrosoguanidine (MNNG)		
2-Nitrofluorene (2-NF)	<u>50</u>	µg/plate TA1538, TA98
9-Aminoacridine (9-AA)	<u>100</u>	µg/plate TA1537

Activation:

2-Anthramine (2-AA)	<u>10</u>	µg/plate all strains
---------------------	-----------	----------------------

### 3. Activation: S9 derived from adult male

<u>x</u> Aroclor 1254	<u>x</u> induced	<u>x</u> rat	<u>x</u> liver
<u>      </u> phenobarbital	<u>      </u> noninduced	<u>      </u> mouse	<u>      </u> lung
<u>      </u> none		<u>      </u> hamster	<u>      </u> other
<u>      </u> other		<u>      </u> other	

<sup>1</sup>Maron, D.M. and Ames, B.N. (1983). Revised Methods for the Salmonella Mutagenicity Test. Mutat. Res. 113:173-215.

# SALMONELLA

Neither the strain of rats nor the source of the S9 were reported.

## S9 mix composition:

<u>Component:</u>	<u>Concentration/mL</u>
Sodium phosphate buffer (pH 7.4)	100 µmoles
Glucose 6-phosphate	5 µmoles
NADP	4 µmoles
KCl	33 µmoles
MgCl <sub>2</sub>	8 µmoles
S9	100 µL

## 4. Test Organism Used: S. typhimurium strains

_____ TA97	<u>x</u>	TA98	<u>x</u>	TA100	_____ TA102	_____ TA104
<u>x</u> TA1535	<u>x</u>	TA1537	<u>x</u>	TA1538		

list any others:

Test organisms were properly maintained: Overnight cultures were prepared from master plates held at 4°C for no more than 1 month. There is some ambiguity about the source of cultures for the master plates. The study author stated that "master plates are prepared from frozen permanents kept at 0°C," but also stated that cultures are stored "in liquid nitrogen at approximately -190°C."

Checked for appropriate genetic markers (rfa mutation, R factor): Yes.

## 5. Test Compound Concentrations Used:

- (a) Preliminary cytotoxicity assay: Up to 10,000 µg/plate; the concentrations tested were not specified.
- (b) Mutation assay: Five doses (50, 150, 500, 1500, and 5,000 µg/plate) were evaluated in triplicate in the presence and absence of S9 activation; all tester strains were used.

## B. TEST PERFORMANCE:

1. Type of Salmonella Assay: x Standard plate test  
 \_\_\_\_\_ Pre-incubation (\_\_\_\_) minutes  
 \_\_\_\_\_ "Prival" modification  
 \_\_\_\_\_ Spot test  
 \_\_\_\_\_ Other (describe)

## 2. Methods:

- (a) Preliminary cytotoxicity/mutation assays: The procedure used for the preliminary cytotoxicity assay was not specified.

For the mutation assay, 0.1 mL of a 16±4-hour culture of the appropriate tester strain, and up to 100 µL of the appropriate test

## SALMONELLA

material dose, solvent, or positive controls were added to tubes containing 2-mL volumes of molten top agar. For the S9-activated tests, 0.5 mL of the S9-cofactor mix was added. Tester strains, and test and control solutions were added as described. The contents of the tubes were mixed, poured over Vogel-Bonner minimal medium E plates, and incubated at 37°C for 48±4 hours. Triplicate plates per strain, per dose, per condition were used for the test compound, solvent, and positive controls.

3. Evaluation Criteria: The test material was considered positive for a particular strain and condition if it caused at least a doubling in the number of revertants of strains TA98 or TA100, or a tripling in the number of revertants of strains TA1535, TA1537, or TA1538.
4. Statistical Analysis: The data were not analyzed statistically.
5. Protocol: See Appendix B.

### C. REPORTED RESULTS

1. Preliminary Cytotoxicity Assay: Neither the data nor the details of the preliminary cytotoxicity assay were reported. The study author stated that the test material was evaluated in strains TA100 and TA1537 at concentrations up to 10,000 µg/plate. The study author further stated that cytotoxicity was observed between 1000 and 10,000 µg/plate. Accordingly, the dose range selected for the mutation assay was 50-5,000 µg/plate +/-S9.
2. Mutation Assay: Representative results from the mutation assay with TREC SPF 15 lotion are presented in Table 1. As shown, solvent control colony counts for strain TA98 +S9 and strain TA100 +/-S9 were well beyond the expected range reported by Maron and Ames<sup>2</sup> (TA98:30-50 revertants/plate; TA100: 120-200 revertants/plate) and either approached the acceptability limit for TA98 (15-75) or exceeded the range for TA100 (60-220) established by other investigators.<sup>3</sup> The results with the nonactivated positive controls (10 µg/plate MNNG and 50 µg/plate 2-NF) further indicated that the higher-than-expected background counts for TA98 and TA100 coincided with a reduction in sensitivity.

Table 1 also shows that there was a marked difference in the number of revertant colonies of TA1538 and TA98 induced by 50 µg/plate 2-NF. Exposure of these strains to comparable doses of 2-NF generally results in the production of approximately equivalent numbers of revertant colonies. While the response (i.e., fold increase) varies as a function of spontaneous mutant counts, the actual number of mutant colonies induced by 2-NF is consistently close. Similarly, marked differences

---

<sup>2</sup>Maron, D.M. and Ames, B.N. (1983).

<sup>3</sup>deSerres, F.J. and Shelby, M.D. 1979. "Recommendations on Data Production and Analysis Using the Salmonella/Microsome Mutagenicity Assay." Environmental Mutagenesis 1:87-92.

## SALMONELLA

between the number of mutant colonies induced in strains TA1535 and TA100 by 10 µg/plate MNNG were observed in this study, even though a similar phenomenon exists for these strains.

Excessive levels of 2-AA (10 µg/plate) were used for the S9-activated positive control. Such high levels can be expected to be cytotoxic, resulting in reduced revertant levels. Finally, the reporting of colony counts >3000/plate is an unacceptable practice. The generally accepted rule, whether counting manually or using an automatic colony counter, is that counts ≥3000 colonies/plate can not be accurately determined.

Based on the above considerations, it was concluded that the study is invalid.

- D. REVIEWERS' DISCUSSION/CONCLUSIONS: The numerous technical deficiencies identified in the review of the study precludes an evaluation of the relevance of the findings with the test material. The study should be repeated using established procedures for the S. typhimurium/mammalian microsome mutagen-icity assay.<sup>4</sup> Additionally, information about the purity of the test material should be submitted with the repeat study.
- E. QUALITY ASSURANCE MEASURES: Was the test performed under GLP? Yes. (A quality assurance statement was signed but not dated.)
- F. CBI APPENDICES: Appendix A, Introduction and Materials and Methods, CBI pp. 7-9; Appendix B, Protocol, CBI pp. 20-23.
- CORE CLASSIFICATION: Unacceptable. The study does not satisfy the data Guideline requirement (§84-2) for genetic effects Category I, Gene Mutations.

---

<sup>4</sup>Maron, D.M. and Ames, B.N. (1983).

TABLE 1: Representative Results of the Salmonella typhimurium/Mammalian Microsome Mutation Assay with TREO SPF 15 Lotion

Substance	Dose/Plate	S9 Activation	Revertants per Plate of Bacterial Tester Strain <sup>a</sup>				
			TA1535	TA1537	TA1538	TA98	TA100
<u>Solvent Control</u>							
Dimethyl sulfoxide	≤100 µL	-	19±6	16±3	29±4	43±2	275±16
	≤100 µL	+	37±4 <sup>b</sup>	13±2	31±7	74±8	370±6
<u>Positive Controls</u>							
N-methyl-N-nitro- N-nitrosoguanidine	10 µg	-	4019±113	--	--	--	1027±123
9-Aminoacridine	100 µg	-	--	641±405	--	--	--
2-Nitrofluorene	50 µg	-	--	--	3044±154	1544±169	--
2-Anthramine	10 µg	+	469±201	1151±13	3549±44	3753±192	3270±103
<u>Test Material</u>							
TREO SPF 15 lotion	500 µg <sup>c</sup>	-	25±5	15±6	30±7	47±5	275±12
	1500 µg	-	18±9	7±2	19±7	36±4	309±46
	5000 µg	-	12±3	0±0	10±5	33±3	286±28
	150 µg <sup>c</sup>	+	28±1	24±3	24±6	109±4	370±23
	500 µg	+	26±11	20±6	28±3	96±13	328±16
	1500 µg	+	28±1	20±6	30±3	90±13	343±32
	5000 µg	+	22±5	20±2	32±10	75±3	306±18

<sup>a</sup>Means and standard deviations of the counts from triplicate plates.

<sup>b</sup>Calculated by our reviewers; erroneously reported by the study author as 27±21.

<sup>c</sup>Results for lower doses (50 and 150 µg/plate -S9, and 50 µg/plate +S9) did not suggest a mutagenic effect.

SALMONELLA

APPENDIX A

INTRODUCTION AND MATERIALS AND METHODS  
CBI pp. 7-9



---

OIL OF CITRONELLA

021901

---

Page \_\_\_\_\_ is not included in this copy.

Pages 9 through 16 are not included.

---

The material not included contains the following type of information:

- ☐ Identity of product inert ingredients.
  - ☐ Identity of product impurities.
  - ☐ Description of the product manufacturing process.
  - ☐ Description of quality control procedures.
  - ☐ Identity of the source of product ingredients.
  - ☐ Sales or other commercial/financial information.
  - ☐ A draft product label.
  - ☐ The product confidential statement of formula.
  - ☐ Information about a pending registration action.
  - ☒ FIFRA registration data.
  - ☐ The document is a duplicate of page(s) \_\_\_\_\_.
  - ☐ The document is not responsive to the request.
- 

The information not included is generally considered confidential by product registrants. If you have any questions, please contact the individual who prepared the response to your request.

---