



325305

**March 4, 2009****SDG Narrative**

**SDG No.: SDG E3YC9**  
**Contract No: EP08W001564**  
**Case No: EPA05-2008-01**

Ms. Terese Van Donsel  
US EPA Region 5  
77 W. Jackson Blvd.  
Chicago, IL 60604

Dear Ms. Donsel,

Enclosed are the results for SDG No.: **SDG E3YC9**. This SDG is associated with Frontier Analytical Laboratory project **5340**. The one soil sample received on 2/11/2009 included the first and last sample of SDG No.: SDG E3YC9. As per your chain of custody instructions, we extracted and analyzed a duplicate sample from the same sample jar. Both soil samples in SDG E3YC9 were extracted and analyzed by EPA Method DLM02.0 for tetra through octa chlorinated dibenzo dioxins and furans.

Please note that this sample was received at a temperature of eighteen degrees Celsius. EPA Method DLM02.0 requires us to notify the SMO if a sample is received above a temperature of ten degrees Celsius. We reported the sample receipt temperature in the sample receipt confirmation paperwork e-mailed to you on 2/12/2009.

In addition, please note that there were matrix interferences throughout several of the dioxin channels; most likely due to high levels of PCBs present in the samples. We were forewarned of the potential high levels of compounds and therefore decided to extract and analyze one gram of each sample. Even with this reduced sample size, the interferences were present. In the case of sample 5340-001-0001-SA (EPA Sample ID: E3YC9), the following is attributed to these interferences and should be noted accordingly. 1,2,3,7,8-PeCDD and 1,2,3,4,7,8-HxCDD have an higher than anticipated EDLs of 12.9 ng/kg and 2.81 ng/kg respectively. 1,2,3,6,7,8 HxCDD has a theoretical ion abundance ratio that is outside the recommended range. Areas from both m/z channels were added and an EMPC is reported. The ion ratio has been qualified with an “\*”, the concentration qualified with the “X” for both the individual target analyte as well as the homologue. In both cases the result has been classified as an EMPC and not as a concentration. In addition, the EMPC value has not been included in the calculation of the TEQ. In the case of sample 5340-002-0001-SA (EPA Sample ID: E3YD0) the following is attributed to these interferences and should be noted accordingly. 2,3,7,8-TCDD, 1,2,3,7,8-PeCDD and 1,2,3,4,7,8-HxCDD have an higher than anticipated EDLs of 1.58 ng/kg, 6.31 ng/kg, and 3.85 ng/kg respectively. 1,2,3,6,7,8 HxCDD as well as some of the PeCDD totals have a theoretical ion abundance ratio that is outside the



recommended range. Areas from both m/z channels were added and an EMPC is reported. The ion ratio has been qualified with an “\*”, the concentration qualified with the “X” for both the individual target analyte as well as the homologue. In both cases the result has been classified as an EMPC and not as a concentration. In addition, the EMPC value has not been included in the calculation of the TEQ.

This SDG used the following GC column; a J&W Scientific/Agilent DB-5 column; 0.25mm ID, 60 meters in length, a Phenyl Arylene polymer coating virtually equivalent to a (5%-Phenyl)-methylpolysiloxane, and a film thickness of 0.25 um. In addition, a second column confirmation was utilized for all samples that had a positive detection of 2,3,7,8-TCDF above the ML. The confirmation column was a J & W Scientific/Agilent DB-225 column; 0.25mm ID, 30 meters in length, a (50%-Cyanopropylphenyl)-dimethylpolysiloxane with a film thickness of 0.25um.

The following are instructions on how to perform the calculations associated with this SDG.

#### ***Response Factor Calculations***

$$RF = \frac{(A_{1r} + A_{2r}) C_r}{(A_{1r} + A_{2r}) C_i}$$

Where:

$C_r$  = The concentration of the labeled compound in the recovery standard.

$C_i$  = The concentration of the labeled compound in the internal standard.

$A_{1r} + A_{2r}$  = The areas of the primary and secondary m/z's for the labeled compound in the recovery standard.

$A_{1i} + A_{2i}$  = The areas of the primary and secondary m/z's for the labeled compound in the internal standard.

RF = The corresponding average response factor derived from the ICAL.

#### ***Concentration Calculations***

$$C_n = \frac{(A_{1n} + A_{2n}) C_i}{(A_{1i} + A_{2i}) RRF(V)}$$

Where:

$C_n$  = The concentration of each PCDD/F in the sample.

$C_i$  = The amount of the labeled compound (internal standard) spiked into the sample.

$A_{1n} + A_{2n}$  = The areas of the primary and secondary m/z's for the PCDD/F.

$A_{1i} + A_{2i}$  = The areas of the primary and secondary m/z's for the labeled compound.

RRF = The corresponding average relative response factor derived from the ICAL.

V = Volume, weight, or amount of sample extracted.



### ***Detection Limit Calculations***

$$DL = \frac{2.5 (H_{n1} + H_{n2}) C_l}{H_{i1} + H_{i2} (RRF) V}$$

Where:

DL = Detection limit of each PCDD/F in the sample.  
H<sub>n1</sub> = Peak to peak noise height of the primary m/z for the PCDD/F.  
H<sub>n2</sub> = Peak to peak noise height of the primary m/z for the PCDD/F.  
C<sub>l</sub> = The amount of the labeled compound (internal standard) spiked into the sample.  
H<sub>i1</sub> = The internal standard peak height of the primary m/z for the PCDD/F.  
H<sub>i2</sub> = The internal standard peak height of the primary m/z for the PCDD/F.  
RRF = The corresponding average relative response factor derived from the ICAL.  
V = Volume, weight, or amount of sample extracted.

### ***Standard Recovery Calculations***

$$\% \text{ Recovery} = \frac{(A_{1r} + A_{2r}) C_r}{(A_{1r} + A_{2r}) (C_l) RF} \times 100$$

Where:

C<sub>r</sub> = The amount of the labeled compound (recovery standard) spiked into the extract.  
C<sub>l</sub> = The amount of the labeled compound (internal standard) spiked into the sample.  
A<sub>1r</sub> + A<sub>2r</sub> = The areas of the primary and secondary m/z's for the labeled compound in the recovery standard.  
A<sub>1l</sub> + A<sub>2l</sub> = The areas of the primary and secondary m/z's for the labeled compound in the internal standard.  
RF = The corresponding average response factor derived from the ICAL.

### ***Percent Recovery***

$$\% \text{ Recovery} = \frac{A_f - A_b}{A_s} \times 100$$

Where:

A<sub>f</sub> = Amount of native found.  
A<sub>b</sub> = Amount of background native.  
A<sub>s</sub> = Amount of native spiked.



*Relative Percent Difference*

$$RPD = \frac{\frac{(\%R_1 - \%R_2)}{(\%R_1 + \%R_2)}}{2} \times 100$$

Where:

$R_1$  = Recovery of analyte in sample one.

$R_2$  = Recovery of analyte in sample two.

If you have any questions regarding SDG No.: SDG E3YC9, Frontier Analytical Laboratory project **5340**, please contact me at (916) 934-0900. The enclosed results are specifically for the samples referenced in this report only. These results meet all NELAC requirements and shall not be reproduced except in full. Thank you for choosing Frontier Analytical Laboratory for your analytical testing needs.

I certify that this data package is in compliance with the terms and conditions of the contract, both technically and for completeness, for other than the conditions detailed above. Release of the data contained in this hardcopy data package has been authorized by the Laboratory Manager or his/her designee, as verified by the following signature.

Sincerely,

A handwritten signature in black ink, appearing to read "Bradley B. Silverbush".

Bradley B. Silverbush  
Director of Operations  
March 4, 2009

1DFA - FORM I-HR CDD-1  
 CDD/CDF SAMPLE DATA SUMMARY  
 HIGH RESOLUTION

SAMPLE No.  
 E3YC9

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564  
 LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: SDG NO.: SDG E3YC9  
 MATRIX: (SOIL/WATER/ASH/TISSUE/OIL) SOIL LAB SAMPLE ID: 5340-001-0001-SA  
 SAMPLE wt/vol: 1.000 (g/mL): g LAB FILE ID: 24FEB09M Sam: 14  
 WATER SAMPLE PREP: SPE (SEPF/SPE)  
 DATE RECEIVED: 11-FEB-09  
 CONCENTRATED EXTRACT VOLUME: 20 (uL)  
 DATE EXTRACTED: 18-FEB-09  
 INJECTION VOLUME: 2 (uL) % SOLIDS/LIPIDS: 75.15  
 DATE ANALYZED: 24-FEB-09  
 GC COLUMN: DB5 ID: 0.25 (mm)  
 DILUTION FACTOR: 1  
 CONCENTRATION UNITS: (pg/L or ng/kg) ng/kg

TARGET ANALYTE	SELECTED IONS	PEAK RT	ION RATIO #	CONCENTRATION	Q	EMPC/EDL
2,3,7,8-TCDD	320/322	27:33	0.67	2.01	J	*
2,3,7,8-TCDF	304/306	26:49	0.68	55.7		*
1,2,3,7,8-PeCDF	340/342	31:40	1.63	155		*
1,2,3,7,8-PeCDD	356/358	NotFnd	*	*	U	12.9
2,3,4,7,8-PeCDF	340/342	32:59	1.54	33.6		*
1,2,3,4,7,8-HxCDF	374/376	37:20	1.26	336		*
1,2,3,6,7,8-HxCDF	374/376	37:32	1.27	154		*
1,2,3,4,7,8-HxCDD	390/392	NotFnd	*	*	U	2.81
1,2,3,6,7,8-HxCDD	390/392	38:53	0.11 *	*	X	116
1,2,3,7,8,9-HxCDD	390/392	39:21	1.07	10.5	J	*
2,3,4,6,7,8-HxCDF	374/376	38:29	1.27	26.1		*
1,2,3,7,8,9-HxCDF	374/376	39:58	1.23	48.0		*
1,2,3,4,6,7,8-HpCDF	408/410	42:24	1.02	646		*
1,2,3,4,6,7,8-HpCDD	424/426	44:18	1.06	132		*
1,2,3,4,7,8,9-HpCDF	408/410	45:14	0.99	184		*
OCDD	458/460	49:54	0.89	10600		*
OCDF	442/444	50:16	0.88	9420		*

NOTE: Concentrations, Estimated Maximum Possible Concentrations (EMPCs), and Estimated Detection Levels (EDLs) for solid samples are calculated on a dry weight basis (except tissues, which are reported on a wet weight basis with % Lipids).

LABELED COMPOUNDS	SELECTED IONS	PEAK RT	ION RATIO #	ION RATIO LIMITS	RECOVERY	
					% REC #	LIMITS
13C-2,3,7,8-TCDD	332/334	27:32	0.77	0.65-0.89	101	25-164
13C-1,2,3,7,8-PeCDD	368/370	33:21	1.75	1.32-1.78	111	25-181
13C-1,2,3,4,7,8-HxCDD	402/404	38:43	1.29	1.05-1.43	86.8	32-141
13C-1,2,3,6,7,8-HxCDD	402/404	38:53	1.29	1.05-1.43	85.9	28-130
13C-1,2,3,4,6,7,8-HpCDD	436/438	44:17	1.06	0.88-1.20	87.7	23-140
13C-OCDD	470/472	49:52	0.92	0.76-1.02	76.7	17-157
13C-2,3,7,8-TCDF	316/318	26:48	0.85	0.65-0.89	101	24-169
13C-1,2,3,7,8-PeCDF	352/354	31:39	1.66	1.32-1.78	101	24-185
13C-2,3,4,7,8-PeCDF	352/354	32:57	1.66	1.32-1.78	99.9	21-178
13C-1,2,3,4,7,8-HxCDF	384/386	37:19	0.53	0.43-0.59	80.2	26-152
13C-1,2,3,6,7,8-HxCDF	384/386	37:31	0.53	0.43-0.59	79.6	26-123
13C-1,2,3,7,8,9-HxCDF	384/386	39:53	0.54	0.43-0.59	93.7	29-147
13C-2,3,4,6,7,8-HxCDF	384/386	38:27	0.53	0.43-0.59	81.6	28-136
13C-1,2,3,4,6,7,8-HpCDF	418/420	42:23	0.44	0.37-0.51	75.3	28-143
13C-1,2,3,4,7,8,9-HpCDF	418/420	45:12	0.44	0.37-0.51	95.5	26-138
13C-OCDF	454/456	50:15	0.91	0.76-1.02	75.4	17-157
37Cl-2,3,7,8-TCDD	328/NA	27:35	NA	NA	114	35-197

# Column to be used to flag values outside (QC) limits.

ANALYST: SC

DATE: 3/21/09

1DFA - FORM I-HR CDD-1  
 CDD/CDF SAMPLE DATA SUMMARY  
 HIGH RESOLUTION

SAMPLE No.  
 E3YD0

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564  
 LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: NA SDG NO.: SDG E3YC9  
 MATRIX: (SOIL/WATER/ASH/TISSUE/OIL) SOIL LAB SAMPLE ID: 5340-002-0001-SA  
 SAMPLE wt/vol: 1.010 (g/mL): g LAB FILE ID: 24FEB09M Sam: 15  
 WATER SAMPLE PREP: SPE (SEPF/SPE) DATE RECEIVED: 11-FEB-09  
 CONCENTRATED EXTRACT VOLUME: 20 (uL) DATE EXTRACTED: 18-FEB-09  
 INJECTION VOLUME: 2 (uL) % SOLIDS/LIPIDS: 75.15 DATE ANALYZED: 24-FEB-09  
 GC COLUMN: DB5 ID: 0.25 (mm)  
 DILUTION FACTOR: 1  
 CONCENTRATION UNITS: (pg/L or ng/kg) ng/kg

TARGET ANALYTE	SELECTED IONS	PEAK RT	ION RATIO #	CONCENTRATION	Q	EMPC/EDL
2,3,7,8-TCDD	320/322	27:34	0.66	1.99	J	*
2,3,7,8-TCDF	304/306	26:48	0.69	39.7		*
1,2,3,7,8-PeCDF	340/342	31:39	1.53	113		*
1,2,3,7,8-PeCDD	356/358	Not Fnd	*	*	U	6.31
2,3,4,7,8-PeCDF	340/342	32:57	1.59	25.7		*
1,2,3,4,7,8-HxCDF	374/376	37:20	1.27	247		*
1,2,3,6,7,8-HxCDF	374/376	37:31	1.27	110		*
1,2,3,4,7,8-HxCDD	390/392	Not Fnd	*	*	U	3.85
1,2,3,6,7,8-HxCDD	390/392	38:52	0.08	*	X	109
1,2,3,7,8,9-HxCDD	390/392	39:18	1.06	8.91	J	*
2,3,4,6,7,8-HxCDF	374/376	38:28	1.23	17.6	J	*
1,2,3,7,8,9-HxCDF	374/376	39:57	1.26	34.1		*
1,2,3,4,6,7,8-HpCDF	408/410	42:24	1.03	475		*
1,2,3,4,6,7,8-HpCDD	424/426	44:18	1.09	99.8		*
1,2,3,4,7,8,9-HpCDF	408/410	45:13	1.04	133		*
OCDD	458/460	49:53	0.89	8050		*
OCDF	442/444	50:16	0.89	6930		*

NOTE: Concentrations, Estimated Maximum Possible Concentrations (EMPCs), and Estimated Detection Levels (EDLs) for solid samples are calculated on a dry weight basis (except tissues, which are reported on a wet weight basis with % Lipids).

Labeled Compounds	SELECTED IONS	PEAK RT	ION RATIO #	ION RATIO LIMITS	RECOVERY	
					% REC #	LIMITS
13C-2,3,7,8-TCDD	332/334	27:32	0.78	0.65-0.89	102	25-164
13C-1,2,3,7,8-PeCDD	368/370	33:20	1.77	1.32-1.78	109	25-181
13C-1,2,3,4,7,8-HxCDD	402/404	38:41	1.30	1.05-1.43	88.1	32-141
13C-1,2,3,6,7,8-HxCDD	402/404	38:51	1.29	1.05-1.43	86.2	28-130
13C-1,2,3,4,6,7,8-HpCDD	436/438	44:17	1.08	0.88-1.20	87.8	23-140
13C-OCDD	470/472	49:52	0.92	0.76-1.02	74.1	17-157
13C-2,3,7,8-TCDF	316/318	26:47	0.84	0.65-0.89	101	24-169
13C-1,2,3,7,8-PeCDF	352/354	31:37	1.67	1.32-1.78	105	24-185
13C-2,3,4,7,8-PeCDF	352/354	32:56	1.66	1.32-1.78	99.4	21-178
13C-1,2,3,4,7,8-HxCDF	384/386	37:17	0.53	0.43-0.59	78.3	26-152
13C-1,2,3,6,7,8-HxCDF	384/386	37:30	0.53	0.43-0.59	77.1	26-123
13C-1,2,3,7,8,9-HxCDF	384/386	39:52	0.54	0.43-0.59	91.1	29-147
13C-2,3,4,6,7,8-HxCDF	384/386	38:26	0.54	0.43-0.59	82.3	28-136
13C-1,2,3,4,6,7,8-HpCDF	418/420	42:23	0.44	0.37-0.51	74.2	28-143
13C-1,2,3,4,7,8,9-HpCDF	418/420	45:12	0.45	0.37-0.51	95.5	26-138
13C-OCDF	454/456	50:15	0.92	0.76-1.02	73.2	17-157
37Cl-2,3,7,8-TCDD	328/NA	27:33	NA	NA	108	35-197

# Column to be used to flag values outside (QC) limits.

ANALYST: 

DATE:  3/6/09

1DFB - FORM I-HR CDD-2  
CDD/CDF TOXICITY EQUIVALENCE SUMMARY  
HIGH RESOLUTION

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: SDG NO.: SDG E3YC9

MATRIX: (SOIL/WATER/ASH/TISSUE/OIL) SOIL LAB SAMPLE ID: 5340-001-0001-SA

SAMPLE wt/vol: 1.000 (g/mL): g LAB FILE ID: 24FEB09M Sam: 14

WATER SAMPLE PREP: SPE (SEPF/SPE) DATE RECEIVED: 11-FEB-09

CONCENTRATED EXTRACT VOLUME: 20 (uL) DATE EXTRACTED: 18-FEB-09

INJECTION VOLUME: 2 (uL) % SOLIDS/LIPIDS: 75.15 DATE ANALYZED: 24-FEB-09

GC COLUMN: DB5 ID: 0.25 (mm) DILUTION FACTOR: 1

CONCENTRATION UNITS: (pg/L or ng/kg) ng/kg

TARGET ANALYTE	CONCENTRATION	TEF*	TEF-ADJUSTED CONCENTRATION
2,3,7,8-TCDD	2.01	x 1.0 =	2.01
2,3,7,8-TCDF	55.7	x 0.1 =	5.57
1,2,3,7,8-PeCDF	155	x 0.05 =	7.75
1,2,3,7,8-PeCDD	*	x 0.5 =	*
2,3,4,7,8-PeCDF	33.6	x 0.5 =	16.8
1,2,3,4,7,8-HxCDF	336	x 0.1 =	33.6
1,2,3,6,7,8-HxCDF	154	x 0.1 =	15.4
1,2,3,4,7,8-HxCDD	*	x 0.1 =	*
1,2,3,6,7,8-HxCDD	*	x 0.1 =	*
1,2,3,7,8,9-HxCDD	10.5	x 0.1 =	1.05
2,3,4,6,7,8-HxCDF	26.1	x 0.1 =	2.61
1,2,3,7,8,9-HxCDF	48.0	x 0.1 =	4.80
1,2,3,4,6,7,8-HpCDF	646	x 0.01 =	6.46
1,2,3,4,6,7,8-HpCDD	132	x 0.01 =	1.32
1,2,3,4,7,8,9-HpCDD	184	x 0.01 =	1.84
OCDD	10600	x 0.001 =	10.6
OCDF	9420	x 0.001 =	9.42

Total = 119

\* TEF - Toxicity Equivalent Factors from EPA/625/3-89/016 March 1989 - Interim Procedures for Estimating Risks Associated with Exposures to Mixtures of Chlorinated Dibenzo-p-Dioxins and Chlorinated Dibenzofurans (CDDs and CDFs) and 1989 Update

ANALYST: [Signature] DATE: 3/2/09

1DFB - FORM I-HR CDD-2  
CDD/CDF TOXICITY EQUIVALENCE SUMMARY  
HIGH RESOLUTION

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: NA SDG NO.: SDG E3YC9

MATRIX: (SOIL/WATER/ASH/TISSUE/OIL) SOIL LAB SAMPLE ID: 5340-002-0001-SA

SAMPLE wt/vol: 1.010 (g/mL): g LAB FILE ID: 24FEB09M Sam: 15

WATER SAMPLE PREP: SPE (SEPF/SPE) DATE RECEIVED: 11-FEB-09

CONCENTRATED EXTRACT VOLUME: 20 (uL) DATE EXTRACTED: 18-FEB-09

INJECTION VOLUME: 2 (uL) % SOLIDS/LIPIDS: 75.15 DATE ANALYZED: 24-FEB-09

GC COLUMN: DB5 ID: 0.25 (mm) DILUTION FACTOR: 1

CONCENTRATION UNITS: (pg/L or ng/kg) ng/kg

TARGET ANALYTE	CONCENTRATION	TEF*	TEF-ADJUSTED CONCENTRATION
2,3,7,8-TCDD	1.99	x 1.0 =	1.99
2,3,7,8-TCDF	39.7	x 0.1 =	3.97
1,2,3,7,8-PeCDF	113	x 0.05 =	5.65
1,2,3,7,8-PeCDD	*	x 0.5 =	*
2,3,4,7,8-PeCDF	25.7	x 0.5 =	12.9
1,2,3,4,7,8-HxCDF	247	x 0.1 =	24.7
1,2,3,6,7,8-HxCDF	110	x 0.1 =	11.0
1,2,3,4,7,8-HxCDD	*	x 0.1 =	*
1,2,3,6,7,8-HxCDD	*	x 0.1 =	*
1,2,3,7,8,9-HxCDD	8.91	x 0.1 =	0.891
2,3,4,6,7,8-HxCDF	17.6	x 0.1 =	1.76
1,2,3,7,8,9-HxCDF	34.1	x 0.1 =	3.41
1,2,3,4,6,7,8-HpCDF	475	x 0.01 =	4.75
1,2,3,4,6,7,8-HpCDD	99.8	x 0.01 =	0.998
1,2,3,4,7,8,9-HpCDF	133	x 0.01 =	1.33
OCDD	8050	x 0.001 =	8.05
OCDF	6930	x 0.001 =	6.93
Total =			88.3

\* TEF - Toxicity Equivalent Factors from EPA/625/3-89/016 March 1989 - Interim Procedures for Estimating Risks Associated with Exposures to Mixtures of Chlorinated Dibenzo-p-Dioxins and Chlorinated Dibenzofurans (CDDs and CDFs) and 1989 Update

ANALYST: E

DATE: 3/2/09

1DFC - FORM I-HR CDD-3  
CDF SECOND COLUMN CONFIRMATION  
HIGH RESOLUTION

SAMPLE No.  
E3YC9

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: SDG NO.: SDG E3YC9

MATRIX: (SOIL/WATER/ASH/TISSUE/OIL) SOIL LAB SAMPLE ID: 5340-001-0001-SA

SAMPLE wt/vol: 1.000 (g/ml) g LAB FILE ID: 04MAR09A Sam: 3

WATER SAMPLE PREP: SPE (SEPF/SPE) DATE RECEIVED: 11-FEB-09

CONCENTRATED EXTRACT VOLUME: 20 (uL) DATE EXTRACTED: 18-FEB-09

INJECTION VOLUME: 2 (uL) % SOLIDS/LIPIDS: 75.15 DATE ANALYZED: 4-MAR-09

GC COLUMN: DB225 ID: 0.25 (mm) DILUTION FACTOR: 1

CONCENTRATION UNITS: (pg/L or ng/kg) ng/kg

ANALYTE	SELECTED IONS	PEAK RT	ION RATIO #	CONCENTRATION	Q	EMPC/EDL
2,3,7,8-TCDF	304/306	19:28	0.77	68.5		6.64

NOTE: Concentrations, Estimated Maximum Possible Concentrations (EMPCs), and Estimated Detection Limits (EDLs) for solid samples are calculated on a dry weight basis (except tissues, which are reported on a wet weight basis with % Lipids).

LABELED COMPOUNDS	SELECTED IONS	PEAK RT	ION RATIO #	ION RATIO LIMITS	% REC #	RECOVERY LIMITS
13C-2,3,7,8-TCDF	316/318	19:26	0.81	0.65-0.89	96.2	24-169

# Column to be used to flag values outside Quality Control (QC) limits.

ANALYST: J

DATE: 3/4/09

1DFC - FORM I-HR CDD-3  
CDF SECOND COLUMN CONFIRMATION  
HIGH RESOLUTION

SAMPLE No.  
E3YD0

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: SDG NO.: SDG E3YC9

MATRIX: (SOIL/WATER/ASH/TISSUE/OIL) SOIL LAB SAMPLE ID: 5340-002-0001-SA

SAMPLE wt/vol: 1.010 (g/ml) g LAB FILE ID: 04MAR09A Sam: 4

WATER SAMPLE PREP: SPE (SEPF/SPE) DATE RECEIVED: 11-FEB-09

CONCENTRATED EXTRACT VOLUME: 20 (uL) DATE EXTRACTED: 18-FEB-09

INJECTION VOLUME: 2 (uL) % SOLIDS/LIPIDS: 75.15 DATE ANALYZED: 4-MAR-09

GC COLUMN: DB225 ID: 0.25 (mm) DILUTION FACTOR: 1

CONCENTRATION UNITS: (pg/L or ng/kg) ng/kg

ANALYTE	SELECTED IONS	PEAK RT	ION RATIO #	CONCENTRATION	Q	EMPC/EDL
2,3,7,8-TCDF	304/306	19:27	0.78	47.7		

NOTE: Concentrations, Estimated Maximum Possible Concentrations (EMPCs), and Estimated Detection Limits (EDLs) for solid samples are calculated on a dry weight basis (except tissues, which are reported on a wet weight basis with % Lipids).



LABELED COMPOUNDS	SELECTED IONS	PEAK RT	ION RATIO #	ION RATIO LIMITS	% REC #	RECOVERY LIMITS
13C-2,3,7,8-TCDF	316/318	19:26	0.80	0.65-0.89	98.1	24-169

# Column to be used to flag values outside Quality Control (QC) limits.

ANALYST: J


DATE: 3/4/09

File:24FEB09M #1-390 Acq:24-FEB-2009 01:35:56 GC El+ Voltage SIR Autospec-Ultima  
319.8965 S:14 BSUB(10000,15,-3.0) PKD(5,5,3,0 10%,100.0 0.00%,F,F) Exp:PCDD  
Sample Text:5340-001-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC




A chromatogram plot with the y-axis labeled 'A' ranging from 0 to 100%. A single, very sharp peak is centered at a retention time of 25.00 minutes. The x-axis has major tick marks at 24.00, 25.00, and 26.00.


File#24FHB09M #1-390 Acq:24-FEB-2009 01:35:56 GC EI+ Voltage SIR Autospec-Ultima  
333.9339 S:14 BSB(10000,15,-3.0) PKD(S,5,3,0,10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:#340-001-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-390 Acq:24-FEB-2009 01:35:56 GC EI+ Voltage SIR Autospec-Ultima  
319.8965 S:14 BSUB(10000,15,-3,0) Exp:PCDD  
Sample Text:5340-001-0001-SA File Text:Frontier Analytical Laboratory



File:24FEB09M #1-424 Acq:24-FEB-2009 01:35:56 GC EI+ Voltage SIR Autospec-Ultima  
355.8546 S:14 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-001-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-424 Acq:24-FEB-2009 01:35:56 GC EI+ Voltage SIR Autospec-Ultima  
357.8517 S:14 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-001-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-424 Acq:24-FEB-2009 01:35:56 GC EI+ Voltage SIR Autospec-Ultima  
367.8949 S:14 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-001-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC





File:24FEB09M #1-424 Acq:24-FEB-2009 01:35:56 GC EI+ Voltage SIR Autospec-Ultima  
369.8919 S:14 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-001-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-424 Acq:24-FEB-2009 01:35:56 GC EI + Voltage SIR Autospec-Ultima  
355.8546 S:14 F:2 BSUB(10000,15,-3,0) Exp:PCDD  
Sample Text:5340-001-0001-SA File Text:Frontier Analytical Laboratory



File:24FEB09M #1-464 Acq:24-FEB-2009 01:35:56 GC EI+ Voltage SIR Autospec-Ultima  
389.8156 S:14 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,FF) Exp:PCDD  
Sample Text:5340-001-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-464 Acq:24-FEB-2009 01:35:56 GC EI+ Voltage SIR Autospec-Ultima  
403.8530 S:14 R:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0.0,0.0%,F,F) EXP:PCDD  
Sample Text:5340-001-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-464 Acq:24-FEB-2009 01:35:56 GC EI+ Voltage SIR Autospec-Ultima  
380.9760 S:14 F:3 Exp:PCDD  
Sample Text:5340-001-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC


File:24FEB09M #1-464 Acq:24-FEB-2009 01:35:56 GC EI+ Voltage SIR Autospec-Ultima  
389.8156 S:14 F:3 BSUB(10000,15,-3.0) Exp:PCDD  
Sample Text:5340-001-0001-SA File Text:Frontier Analytical Laboratory




File:24FEB09M #1-541 Acq:24-FEB-2009 01:35:56 GC EI+ Voltage SIR Autospec-Ultima  
423.7767 S:14 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100,0,0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-001-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-541 Acq:24-FEB-2009 01:35:56 GC EI+ Voltage SIR Autospec-Ultima  
425.7737 S:14 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100,0,0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-001-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-541 Acq:24-FEB-2009 01:35:56 GC EI+ Voltage SIR Autospec-Ultima  
437.8140 S:14 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100,0,0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-001-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-347 Acq:24-FEB-2009 01:35:56 GC EI+ Voltage SIR Autospec-Ultima  
457.7377 S:14 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-001-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-347 Acq:24-FEB-2009 01:35:56 GC EI+ Voltage SIR Autospec-Ultima  
469.7780 S:14 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-001-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-347 Acq:24-FEB-2009 01:35:56 GC EI+ Voltage SIR Autospec-Ultima  
454.9728 S:14 F:5 Exp:PCDD  
Sample Text:5340-001-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-390 Acq:24-FEB-2009 01:35:56 GC EI+ Voltage SIR Autospec-Ultima  
303.9016 S:14 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-001-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-390 Acq:24-FEB-2009 01:35:56 GC EI+ Voltage SIR Autospec-Ultima  
305.8987 S:14 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-001-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-390 Acq:24-FEB-2009 01:35:56 GC EI+ Voltage SIR Autospec-Ultima  
317.9389 S:14 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-001-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC

100 %  
50 %  
0 %




File:24FEB09M #1-390 Acq:24-FEB-2009 01:35:56 GC EI+ Voltage SIR Autospec-Ultima  
375.8364 S:14 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-001-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC






File:24FEB09M #1-390 Acq:24-FEB-2009 01:35:56 GC EI+ Voltage SIR Autospec-Ultima  
303.9016 S:14 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-001-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC





000025 of 000279




File:24FEB09M #1-390 Acq:24-FEB-2009 01:35:56 GC EI+ Voltage SIR Autospec-Ultima  
339.8597 S:14 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-001-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC





File:24FEB09M #1-424 Acq:24-FEB-2009 01:35:56 GC EI+ Voltage SIR Autospec-Ultima  
 339.8597 S:14 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
 Sample Text:5340-001-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC







File:24FEB09M #1-424 Acq:24-FEB-2009 01:35:56 GC EI+ Voltage SIR Autospec-Ultima  
339.8397 S:14 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-001-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC







File:24FEB09M #1-464 Acq:24-FEB-2009 01:35:56 GC EI+ Voltage SIR Autospec-Ultima  
 373.8207 S:14 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
 Sample Text:5340-001-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-541 Acq:24-FEB-2009 01:35:56 GC EI+ Voltage SIR Autospec-Ultima  
407.7818 S:14 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-001-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-347 Acq:24-FEB-2009 01:35:56 GC EI+ Voltage SIR Autospec-Ultima  
441.7428 S:14 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-001-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC




File:04MAR09A #1-1388 Acq: 4-MAR-2009 10:56:46 GC El+ Voltage SIR 70S  
303.9016 S:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:TCDF  
Sample Text:5340-001-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:BS




File:24FEB09M #1-390 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
319.8965 S:15 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD

Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-390 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
321.8936 S:15 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD

Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-390 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
327.8847 S:15 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD

Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-390 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
331.9368 S:15 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD

Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-390 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
319.8965 S:15 BSUB(10000,15,-3,0) Exp:PCDD  
Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory




File:24FEB09M #1-425 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
355.8546 S:15 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-425 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
357.8517 S:15 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-425 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
369.8919 S:15 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-425 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
355.8546 S:15 F:2 BSUB(10000,15,-3.0) Exp:PCDD  
Sample Text:53340-002-0001-SA File Text:Frontier Analytical Laboratory

A6.05E4

1.1E4

9.6E3

8.5E3

7.5E3

6.4E3

5.3E3


4.3E3

3.2E3

2.1E3


1.1E3

0.0E0







File:24FEB09M #1-425 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
357.8517 S:15 F:2 BSUB(10000,15,-3.0) Exp:PCDD  
Sample Text:53340-002-0001-SA File Text:Frontier Analytical Laboratory


A7.13E4



File:24FEB09M #1-464 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
389.8156 S:15 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-464 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
389.8156 S:15 F:3 BSUB(10000,15,-3.0) Exp:PCDD  
Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory  
A3.15E5



File:24FEB09M #1-541 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
423.7767 S:15 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100.0,0.00%,R,F) Exp:PCDD

Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC

42:55 A4.78E5

A4.53E5

44:18  
A.53E5

9.9E4

4.9E4

0.0E0



Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC

42:00 A3.29E5

A5.35E5

43:29  
A8.28E5

A2.59E5

44:18  
A4.17E5


44:54  
A8.32E4

45:21  
A6.07E5

1.8E5

9.2E4

0.0E0



Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC

42:00 A3.29E5

A5.51E4

43:50  
A2.59E5


44:00  
A4.17E5

44:18  
A4.53E5

1.9E6

9.4E5

0.0E0



Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC

42:00 A3.29E5

A5.51E4

43:50  
A2.59E5

44:00  
A4.17E5

44:18  
A4.53E5

1.7E6

8.6E5

0.0E0



Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC

42:00 A3.29E5

A5.51E4

43:50  
A2.59E5


44:00  
A4.17E5

44:18  
A4.53E5

1.7E6

8.6E5

0.0E0



Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC

42:00 A3.29E5

A5.51E4

43:50  
A2.59E5

44:00  
A4.17E5

44:18  
A4.53E5


1.7E6

8.6E5

0.0E0


File:24FEB09M #1-347 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
457.7377 S:15 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC

49:53  
A1.89E7




File:24FEB09M #1-347 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
459.7348 S:15 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC

49:53  
A2.12E7




File:24FEB09M #1-347 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
469.7780 S:15 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC

49:52  
A1.01E7



File:24FEB09M #1-347 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
471.7750 S:15 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC

49:52  
A1.09E7




File:24FEB09M #1-347 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
454.9728 S:15 F:5 Exp:PCDD  
Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC

100 %  
50  
0


49:00 49:30 50:04 50:34 50:50 51:09 51:24 51:53 52:21 Time

1.4E7  
6.8E6  
0.0E0


File:24FEB09M #1-390 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
303,9016 S:15 BSUB(10000,15,-3,0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-390 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
305,8987 S:15 BSUB(10000,15,-3,0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-390 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
315,9419 S:15 BSUB(10000,15,-3,0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-390 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
317,9389 S:15 BSUB(10000,15,-3,0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC





File:24FEB09M #1-390 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
375,8364 S:15 BSUB(10000,15,-3,0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-390 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
303.9016 S:15 BSUB(10000,15,-3.0) Exp:PCDD  
Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory




File#24FEB09M #1-390 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
339.8597 S:15 BSB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0 0.00%,FF) Exp:PCDD  
Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-425 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
339.8597 S:15 F:2 BSUB(10000,15,-3,0) PKD(5,5,3,0,10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-425 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
341.8568 S:15 F:2 BSUB(10000,15,-3.0) PKD(5,5,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-425 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
353.8970 S:15 R:2 BSUB(10000,15,-3.0) PKD(5.5,3.0)10%100.0,0.00% FF) Exp:PCDDF







File:24FEB09M #1-425 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
409.7974 S:15 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC





File:24FEB09M #1-425 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SR AutoSpec-Ultima  
339.8597 S:15 F:2 BSUB(10000,15,-3.0) Exp:PCDD  
Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory  
100 %  
A1.42E6



File:24FEB09M #1-464 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
373.8207 S:15 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-541 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
407.7818 S:15 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-347 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
441.7428 S:15 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100,0,0,00%,F,F) Exp:PCDD  
Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC

50:16  
A1.95E7




File:24FEB09M #1-347 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
453.7831 S:15 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100,0,0,00%,F,F) Exp:PCDD  
Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC

50:16  
A2.19E7  
3.5E6  
1.7E6  
0.0E0



File:24FEB09M #1-347 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
455.7801 S:15 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100,0,0,00%,F,F) Exp:PCDD  
Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC




50:14  
A1.61E7  
50:15  
A1.48E7  
2.3E6  
1.2E6  
0.0E0



File:24FEB09M #1-347 Acq:24-FEB-2009 02:31:11 GC EI+ Voltage SIR Autospec-Ultima  
513.6775 S:15 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100,0,0,00%,F,F) Exp:PCDD  
Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:TC

48:25  
A2.17E3  
48:47  
A614.41  
48:39  
A276.33  
48:57  
A622.72  
49:47  
A785.39  
50:22  
A600.17  
50:50  
A751.61  
51:19  
A354.74  
51:50  
A203.79  
52:05  
A703.08  
1.4E3  
6.8E2  
0.0E0

File:04MAR09A #1-1388 Acq: 4-MAR-2009 11:31:50 GC EI+ Voltage SIR 70S  
 313.9016 S:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:TCDF  
 Sample Text:5340-002-0001-SA File Text:Frontier Analytical Laboratory Analyst ID:BS



FAL ID: 5340-001-0001-SA      Filename: 24FEB09M      Sam:14      Acquired: 24-FEB-09 01:35:56      iCal: PCDDFAL3-2-5-09  
 Client ID: E3YC9      ConCal: ST022409M3      EndCal: ST022409M4      EPA  
 Results: 5340      GC Column: DB5      Amount: 1.000      NATO 1989 Tox: 121      WHO 1998 Tox: 103      WHO 2005 Tox: 92.4  
 Name      Resp      RA      RT      RRF      Conc      Qual      Fac Noise-1      Noise-2      DL

2,3,7,8-TCDD	2.97e+04	0.67 y	27:33	1.04	2.01	J	2.50	-	-	*	
1,2,3,7,8-PeCDD	*	* n	NotFnd	0.90	*	U	2.50	2620	8480	12.9	
1,2,3,4,7,8-HxCDD	*	* n	NotFnd	1.43	*	U	2.50	1430	2130	2.81	
1,2,3,6,7,8-HxCDD	2.21e+05	0.11 n	38:53	1.02	116	X	2.50	-	-	*	
1,2,3,7,8,9-HxCDD	1.17e+05	1.07 y	39:21	1.05	10.5	J	2.50	-	-	*	
1,2,3,4,6,7,8-HpCDD	1.18e+06	1.06 y	44:18	0.99	132		2.50	-	-	*	
OCDD	5.61e+07	0.89 y	49:54	0.94	10600		2.50	-	-	*	
2,3,7,8-TCDF	1.74e+06	0.68 y	26:49	1.34	55.7		2.50	-	-	*	
1,2,3,7,8-PeCDF	3.28e+06	1.63 y	31:40	1.03	155		2.50	-	-	*	
2,3,4,7,8-PeCDF	5.64e+05	1.54 y	32:59	0.81	33.6		2.50	-	-	*	
1,2,3,4,7,8-HxCDF	5.66e+06	1.26 y	37:20	1.35	336		2.50	-	-	*	
1,2,3,6,7,8-HxCDF	2.48e+06	1.27 y	37:32	1.25	154		2.50	-	-	*	
2,3,4,6,7,8-HxCDF	4.85e+05	1.27 y	38:29	1.53	26.1		2.50	-	-	*	
1,2,3,7,8,9-HxCDF	6.20e+05	1.23 y	39:58	1.19	48.0		2.50	-	-	*	
1,2,3,4,6,7,8-HpCDF	7.14e+06	1.02 y	42:24	1.10	646		2.50	-	-	*	
1,2,3,4,7,8,9-HpCDF	1.98e+06	0.99 y	45:14	1.02	184		2.50	-	-	*	
OCDF	5.93e+07	0.88 y	50:16	0.76	9420		2.50	-	-	*	
13C-2,3,7,8-TCDD	2.85e+07	0.77 y	27:32	0.94	2010				Rec		
13C-1,2,3,7,8-PeCDD	2.52e+07	1.75 y	33:21	0.75	2210					101	
13C-1,2,3,4,7,8-HxCDD	2.06e+07	1.29 y	38:43	1.03	1740					111	
13C-1,2,3,6,7,8-HxCDD	2.14e+07	1.29 y	38:53	1.08	1720					86.8	
13C-1,2,3,4,6,7,8-HpCDD	1.81e+07	1.06 y	44:17	0.89	1750					85.9	
13C-OCDD	2.25e+07	0.92 y	49:52	0.64	3070					87.7	
13C-2,3,7,8-TCDF	4.66e+07	0.85 y	26:48	0.86	2010					76.7	
13C-1,2,3,7,8-PeCDF	4.11e+07	1.66 y	31:39	0.76	2020					101	
13C-2,3,4,7,8-PeCDF	4.15e+07	1.66 y	32:57	0.78	2000					101	
13C-1,2,3,4,7,8-HxCDF	2.49e+07	0.53 y	37:19	1.34	1600					99.9	
13C-1,2,3,6,7,8-HxCDF	2.58e+07	0.53 y	37:31	1.40	1590					80.2	
13C-2,3,4,6,7,8-HxCDF	2.43e+07	0.53 y	38:27	1.29	1630					79.6	
13C-1,2,3,7,8,9-HxCDF	2.18e+07	0.54 y	39:53	1.01	1870					81.6	
13C-1,2,3,4,6,7,8-HpCDF	2.01e+07	0.44 y	42:23	1.16	1510					93.7	
13C-1,2,3,4,7,8,9-HpCDF	2.11e+07	0.44 y	45:12	0.96	1910					75.3	
13C-OCDF	3.30e+07	0.91 y	50:15	0.95	3010					95.5	
37Cl-2,3,7,8-TCDD	1.22e+07		27:35	0.88	915					75.4	
13C-1,2,3,4-TCDD	3.03e+07	0.79 y	26:58	-	81.4					114	
13C-1,2,3,4-TCDF	5.36e+07	0.86 y	25:43	-	89.5						
13C-1,2,3,7,8,9-HxCDD	2.31e+07	1.28 y	39:20	-	115						
Total Tetra-Dioxins	1.87e+05		24:34	1.04	12.6		2.50	-	-	*	6
Total Penta-Dioxins	2.94e+05		30:24	0.90	25.9		2.50	-	-	*	2
Total Hexa-Dioxins	1.69e+06		36:16	1.16	142	X	2.50	-	-	*	5
Total Hepta-Dioxins	2.53e+06		42:56	0.99	283		2.50	-	-	*	2
Total Tetra-Furans	1.54e+07		23:15	1.34	494		2.50	-	-	*	18
1st Fn. Tot Penta-Furans	1.25e+06		28:38	0.92	65.8		2.50	-	-	*	PeCDF 1
Total Penta-Furans	1.11e+07		30:14	0.92	569		2.50	-	-	*	635 14
Total Hexa-Furans	1.34e+07		35:43	1.34	822		2.50	-	-	*	15
Total Hepta-Furans	1.17e+07		42:24	1.06	1060		2.50	-	-	*	4

Analyst: R

Date: 3/3/09

Totals class: Total Tetra-Dioxins

Entry #: 38

Run: 11 File: 24FEB09M S: 14 I: 1 F: 1  
Acquired: 24-FEB-09 01:35:56

Total Concentration: 12.6 Unnamed Concentration: 10.634

RT	m1 Resp	m2 Resp RA	Resp	Concentration	Name
24:34	1.62e+04	2.29e+04 0.71 y	3.92e+04	2.65	
24:51	9.50e+03	1.42e+04 0.67 y	2.37e+04	1.60	
25:59	1.15e+04	1.38e+04 0.83 y	2.52e+04	1.70	
26:07	1.37e+04	1.83e+04 0.75 y	3.20e+04	2.16	
27:33	1.19e+04	1.78e+04 0.67 y	2.97e+04	2.01	2,3,7,8-TCDD
27:52	1.48e+04	2.26e+04 0.65 y	3.74e+04	2.53	

Totals class: Total Penta-Dioxins              Entry #: 39

Run: 11            File: 24FEB09M              S: 14 I: 1 F: 2  
Acquired: 24-FEB-09 01:35:56

Total Concentration: 25.9              Unnamed Concentration: 25.908

RT	m1 Resp	m2 Resp RA	Resp	Concentration	Name
30:24	9.94e+04	7.47e+04 1.33 y	1.74e+05	15.3	
31:40	7.37e+04	4.66e+04 1.58 y	1.20e+05	10.6	

Totals class: Total Hexa-Dioxins

Entry #: 40

Run: 11 File: 24FEB09M S: 14 I: 1 F: 3  
Acquired: 24-FEB-09 01:35:56

Total Concentration: 254 Unnamed Concentration: 127.000

RT	m1 Resp	m2 Resp RA	Resp	Concentration	Name
36:16	2.03e+05	1.88e+05 1.08 y	3.91e+05	32.0	
37:12	6.80e+04	1.02e+05 0.67 n	1.23e+05	14.0	
37:37	4.66e+05	5.20e+05 0.90 n	8.43e+05	81.0	
38:53	1.22e+05	1.14e+06 0.11 n	2.21e+05	116	1,2,3,6,7,8-HxCDD
39:21	6.03e+04	5.63e+04 1.07 y	1.17e+05	10.5	1,2,3,7,8,9-HxCDD

Totals class: Total Hepta-Dioxins              Entry #: 41

Run: 11            File: 24FEB09M              S: 14 I: 1 F: 4  
Acquired: 24-FEB-09 01:35:56

Total Concentration: 283              Unnamed Concentration: 151.330

RT	m1 Resp	m2 Resp RA	Resp	Concentration	Name
42:56	6.49e+05	7.02e+05 0.92 y	1.35e+06	151	
44:18	6.06e+05	5.71e+05 1.06 y	1.18e+06	132	1,2,3,4,6,7,8-HpCDD

Totals class: Total Tetra-Furans

Entry #: 42

Run: 11 File: 24FEB09M S: 14 I: 1 F: 1  
Acquired: 24-FEB-09 01:35:56

Total Concentration: 494 Unnamed Concentration: 438.762

RT	m1 Resp	m2 Resp RA	Resp	Concentration	Name
23:15	1.83e+04	2.53e+04 0.73 y	4.36e+04	1.40	
23:36	7.16e+04	1.03e+05 0.70 y	1.75e+05	5.60	
23:57	1.12e+05	1.69e+05 0.66 y	2.81e+05	9.03	
24:22	5.83e+05	8.24e+05 0.71 y	1.41e+06	45.1	
24:36	8.84e+05	1.29e+06 0.68 y	2.18e+06	69.8	
24:54	3.42e+05	5.01e+05 0.68 y	8.43e+05	27.0	
25:07	2.20e+05	3.18e+05 0.69 y	5.38e+05	17.3	
25:26	3.64e+05	5.34e+05 0.68 y	8.99e+05	28.8	
25:37	1.64e+05	2.49e+05 0.66 y	4.13e+05	13.2	
25:43	4.99e+05	7.39e+05 0.68 y	1.24e+06	39.7	
26:07	1.32e+06	1.90e+06 0.70 y	3.22e+06	103	
26:19	1.67e+05	2.36e+05 0.71 y	4.03e+05	12.9	
26:27	5.38e+04	7.36e+04 0.73 y	1.27e+05	4.09	
26:44	4.28e+05	6.18e+05 0.69 y	1.05e+06	33.6	
26:49	7.05e+05	1.03e+06 0.68 y	1.74e+06	55.7	2,3,7,8-TCDF
27:10	1.26e+05	1.83e+05 0.69 y	3.09e+05	9.92	
27:23	7.50e+04	1.13e+05 0.67 y	1.88e+05	6.02	
28:43	1.69e+05	2.04e+05 0.83 y	3.73e+05	12.0	

Totals class: 1st Fn. Tot Penta-Furans      Entry #: 43

Run: 11      File: 24FEB09M      S: 14 I: 1 F: 1  
Acquired: 24-FEB-09 01:35:56

Total Concentration: 65.8      Unnamed Concentration: 65.769

RT	m1 Resp	m2 Resp RA	Resp	Concentration	Name
28:38	7.75e+05	4.73e+05 1.64 y	1.25e+06	65.8	

Totals class: Total Penta-Furans

Entry #: 44

Run: 11 File: 24FEB09M S: 14 I: 1 F: 2  
Acquired: 24-FEB-09 01:35:56

Total Concentration: 569 Unnamed Concentration: 380.423

RT	m1 Resp	m2 Resp RA	Resp	Concentration	Name
30:14	2.37e+05	1.62e+05 1.47 y	3.99e+05	21.0	
30:24	2.07e+06	1.28e+06 1.62 y	3.35e+06	177	
30:53	4.24e+04	3.13e+04 1.35 y	7.37e+04	3.88	
31:06	1.29e+05	8.78e+04 1.47 y	2.17e+05	11.4	
31:14	2.66e+05	1.67e+05 1.60 y	4.32e+05	22.8	
31:24	3.72e+05	2.29e+05 1.63 y	6.01e+05	31.7	
31:40	2.04e+06	1.25e+06 1.63 y	3.28e+06	155	1,2,3,7,8-PeCDF
32:02	7.46e+05	4.81e+05 1.55 y	1.23e+06	64.7	
32:42	3.22e+04	2.32e+04 1.39 y	5.54e+04	2.92	
32:48	6.29e+04	3.61e+04 1.74 y	9.90e+04	5.22	
32:59	3.42e+05	2.22e+05 1.54 y	5.64e+05	33.6	2,3,4,7,8-PeCDF
33:01	3.25e+05	1.93e+05 1.68 y	5.18e+05	27.3	
33:54	5.64e+04	3.27e+04 1.73 y	8.91e+04	4.70	
34:22	9.30e+04	6.32e+04 1.47 y	1.56e+05	8.23	

Totals class: Total Hexa-Furans

Entry #: 45

Run: 11 File: 24FEB09M S: 14 I: 1 F: 3  
Acquired: 24-FEB-09 01:35:56

Total Concentration: 822 Unnamed Concentration: 257.688

RT	m1 Resp.	m2 Resp RA	Resp	Concentration	Name
35:43	2.19e+05	1.81e+05 1.21 y	4.00e+05	24.8	
35:49	1.09e+06	8.54e+05 1.27 y	1.94e+06	120	
35:59	2.07e+04	1.80e+04 1.15 y	3.87e+04	2.39	
36:05	2.70e+04	1.97e+04 1.37 y	4.67e+04	2.89	
36:17	6.94e+04	5.28e+04 1.32 y	1.22e+05	7.55	
36:34	5.01e+05	3.97e+05 1.26 y	8.97e+05	55.5	
36:53	2.55e+04	1.93e+04 1.32 y	4.48e+04	2.77	
37:10	1.36e+05	1.05e+05 1.29 y	2.41e+05	14.9	
37:20	3.16e+06	2.51e+06 1.26 y	5.66e+06	336	1,2,3,4,7,8-HxCDF
37:32	1.39e+06	1.09e+06 1.27 y	2.48e+06	154	1,2,3,6,7,8-HxCDF
37:47	8.37e+04	6.69e+04 1.25 y	1.51e+05	9.31	
37:58	4.66e+04	3.74e+04 1.25 y	8.40e+04	5.20	
38:10	1.09e+05	8.96e+04 1.22 y	1.99e+05	12.3	
38:29	2.71e+05	2.14e+05 1.27 y	4.85e+05	26.1	2,3,4,6,7,8-HxCDF
39:58	3.42e+05	2.79e+05 1.23 y	6.20e+05	48.0	1,2,3,7,8,9-HxCDF

Totals class: Total Hepta-Furans

Entry #: 46

Run: 11 File: 24FEB09M S: 14 I: 1 F: 4  
Acquired: 24-FEB-09 01:35:56

Total Concentration: 1060 Unnamed Concentration: 234.757

RT	m1 Resp	m2 Resp RA	Resp	Concentration	Name
42:24	3.61e+06	3.53e+06 1.02 y	7.14e+06	646	1,2,3,4,6,7,8-HpCDF
42:57	2.52e+05	2.18e+05 1.16 y	4.70e+05	42.8	
43:14	1.07e+06	1.04e+06 1.02 y	2.11e+06	192	
45:14	9.89e+05	9.94e+05 0.99 y	1.98e+06	184	1,2,3,4,7,8,9-HpCDF

FAL ID: 5340-001-0001-SA      Filename: 04MAR09A      Sam:3      Acquired: 4-MAR-09 10:56:46      iCal: TCDFAL1-1-13-09  
Client ID: E3YC9      ConCal: ST030409A1 EndCal: ST030409A2  
Results:      GC Column: DB225      Amount: 1.000

Name	Resp	RA	RT	RRF	Conc	Qual	Fac	Noise	DL	#Hom	Rec
2,3,7,8-TCDF	2.39e+06	0.77	y	19:28	1.03	68.5	2.50	-	-	1	96.2
13C-2,3,7,8-TCDF	6.81e+07	0.81	y	19:26	0.87	1920					
13C-1,2,3,4-TCDF	8.13e+07	0.81	y	16:53	-	43.7					

Analyst: J Date: 3/4/09

FAL ID: 5340-002-0001-SA      Filename: 24FEB09M      Sam:15      Acquired: 24-FEB-09 02:31:11      iCal: PCDDFAL3-2-5-09      EPA  
 Client ID: E3YDO      ConCal: ST022409M3      EndCal: ST022409M4      *TP*  
 Results: 5340      GC Column: DB5      Amount: 1.010      NATO 1989 Tox: 89.8      WHO 1998 Tox: 76.3      WHO 2005 Tox: 71.9  
 Name      Resp      RA      RT      RRF      Conc      Qual      Fac Noise-1      Noise-2      DL  
 2,3,7,8-TCDD      2.92e+04      0.66 y      27:34      1.04      1.99      J      2.50      -      -      \*  
 1,2,3,7,8-PeCDD      \*      \* n NotFnd      0.90      \*      U      2.50      1390      4330      6.31  
 1,2,3,4,7,8-HxCDD      \*      \* n NotFnd      1.43      \*      U      2.50      776      3860      3.85  
 1,2,3,6,7,8-HxCDD      1.61e+05      0.08 n      38:52      1.02      109      X      2.50      -      -      \*  
 1,2,3,7,8,9-HxCDD      9.68e+04      1.06 y      39:18      1.05      8.91      J      2.50      -      -      \*  
 1,2,3,4,6,7,8-HpCDD      8.70e+05      1.09 y      44:18      0.99      99.8      2.50      -      -      \*  
 OCDD      4.01e+07      0.89 y      49:53      0.94      8050      2.50      -      -      \*  
 2,3,7,8-TCDF      1.21e+06      0.69 y      26:48      1.34      39.7      2.50      -      -      \*  
 1,2,3,7,8-PeCDF      2.41e+06      1.53 y      31:39      1.03      113      2.50      -      -      \*  
 2,3,4,7,8-PeCDF      4.18e+05      1.59 y      32:57      0.81      25.7      2.50      -      -      \*  
 1,2,3,4,7,8-HxCDF      3.96e+06      1.27 y      37:20      1.35      247      2.50      -      -      \*  
 1,2,3,6,7,8-HxCDF      1.68e+06      1.27 y      37:31      1.25      110      2.50      -      -      \*  
 2,3,4,6,7,8-HxCDF      3.22e+05      1.23 y      38:28      1.53      17.6      J      2.50      -      -      \*  
 1,2,3,7,8,9-HxCDF      4.18e+05      1.26 y      39:57      1.19      34.1      2.50      -      -      \*  
 1,2,3,4,6,7,8-HpCDF      5.04e+06      1.03 y      42:24      1.10      475      2.50      -      -      \*  
 1,2,3,4,7,8,9-HpCDF      1.40e+06      1.04 y      45:13      1.02      133      2.50      -      -      \*  
 OCDF      4.13e+07      0.89 y      50:16      0.76      6930      2.50      -      -      \*  
 13C-2,3,7,8-TCDD      2.78e+07      0.78 y      27:32      0.94      2010      Rec      102  
 13C-1,2,3,7,8-PeCDD      2.40e+07      1.77 y      33:20      0.75      2150      109  
 13C-1,2,3,4,7,8-HxCDD      2.02e+07      1.30 y      38:41      1.03      1740      88.1  
 13C-1,2,3,6,7,8-HxCDD      2.07e+07      1.29 y      38:51      1.08      1710      86.2  
 13C-1,2,3,4,6,7,8-HpCDD      1.75e+07      1.08 y      44:17      0.89      1740      87.8  
 13C-OCDD      2.10e+07      0.92 y      49:52      0.64      2940      74.1  
 13C-2,3,7,8-TCDF      4.50e+07      0.84 y      26:47      0.86      2000      101  
 13C-1,2,3,7,8-PeCDF      4.09e+07      1.67 y      31:37      0.76      2070      105  
 13C-2,3,4,7,8-PeCDF      3.98e+07      1.66 y      32:56      0.78      1970      99.4  
 13C-1,2,3,4,7,8-HxCDF      2.35e+07      0.53 y      37:17      1.34      1550      78.3  
 13C-1,2,3,6,7,8-HxCDF      2.41e+07      0.53 y      37:30      1.40      1530      77.1  
 13C-2,3,4,6,7,8-HxCDF      2.36e+07      0.54 y      38:26      1.29      1630      82.3  
 13C-1,2,3,7,8,9-HxCDF      2.04e+07      0.54 y      39:52      1.01      1800      91.1  
 13C-1,2,3,4,6,7,8-HpCDF      1.92e+07      0.44 y      42:23      1.16      1470      74.2  
 13C-1,2,3,4,7,8,9-HpCDF      2.04e+07      0.45 y      45:12      0.96      1890      95.5  
 13C-OCDF      3.09e+07      0.92 y      50:15      0.95      2900      73.2  
 37Cl-2,3,7,8-TCDD      1.11e+07      27:33      0.88      853      108  
 13C-1,2,3,4-TCDD      2.93e+07      0.79 y      26:57      -      78.0  
 13C-1,2,3,4-TCDF      5.16e+07      0.85 y      25:42      -      85.3  
 13C-1,2,3,7,8,9-HxCDD      2.23e+07      1.31 y      39:18      -      110  
 Total Tetra-Dioxins      1.11e+05      24:32      1.04      7.61      2.50      -      -      \*      4  
 Total Penta-Dioxins      9.95e+04      30:24      0.90      9.11      X      2.50      -      -      \*      1  
 Total Hexa-Dioxins      1.17e+06      36:15      1.16      99.8      X      2.50      -      -      \*      5  
 Total Hepta-Dioxins      1.88e+06      42:55      0.99      216      2.50      -      -      \*      2  
 Total Tetra-Furans      1.07e+07      23:12      1.34      351      2.50      -      -      \*      19  
 1st Fn. Tot Penta-Furans      6.10e+05      28:35      0.92      32.6      2.50      -      -      \* PeCDF      1  
 Total Penta-Furans      7.69e+06      30:14      0.92      398      2.50      -      -      \* 431      11  
 Total Hexa-Furans      8.85e+06      35:44      1.34      569      2.50      -      -      \*      13  
 Total Hepta-Furans      8.25e+06      42:24      1.06      777      2.50      -      -      \*      4

Analyst: *[Signature]*

Date: *3/3/09*

Totals class: Total Tetra-Dioxins

Entry #: 38

Run: 12 File: 24FEB09M S: 15 I: 1 F: 1  
Acquired: 24-FEB-09 02:31:11

Total Concentration: 7.61 Unnamed Concentration: 5.613

RT	m1 Resp	m2 Resp RA	Resp	Concentration	Name
24:32	1.41e+04	1.66e+04 0.85 y	3.07e+04	2.10	
24:49	9.70e+03	1.12e+04 0.86 y	2.09e+04	1.43	
27:34	1.16e+04	1.76e+04 0.66 y	2.92e+04	1.99	2,3,7,8-TCDD
27:52	1.26e+04	1.79e+04 0.71 y	3.05e+04	2.08	

Totals class: Total Penta-Dioxins

Entry #: 39

Run: 12 File: 24FEB09M

S: 15 I: 1 F: 2

Acquired: 24-FEB-09 02:31:11

Total Concentration: 12.1

Unnamed Concentration: 12.100

RT	m1 Resp	m2 Resp RA	Resp	Concentration	Name
30:24	6.05e+04	7.13e+04	0.85 n	1.32e+05	12.1

Totals class: Total Hexa-Dioxins

Entry #: 40

Run: 12 File: 24FEB09M S: 15 I: 1 F: 3  
Acquired: 24-FEB-09 02:31:11

Total Concentration: 209 Unnamed Concentration: 91.000

RT	m1 Resp	m2 Resp RA	Resp	Concentration	Name
36:15	1.25e+05	1.16e+05 1.07 y	2.41e+05	20.0	
37:11	5.63e+04	7.77e+04 0.72 n	1.34e+05	11.2	
37:36	3.15e+05	4.00e+05 0.79 n	7.15e+05	59.8	
38:52	8.90e+04	1.07e+06 0.08 n	1.16e+06	109	1,2,3,6,7,8-HxCDD
39:18	4.99e+04	4.69e+04 1.06 y	9.68e+04	8.91	1,2,3,7,8,9-HxCDD

Totals class: Total Hepta-Dioxins

Entry #: 41

Run: 12 File: 24FEB09M S: 15 I: 1 F: 4  
Acquired: 24-FEB-09 02:31:11

Total Concentration: 216 Unnamed Concentration: 116.116

RT	m1 Resp	m2 Resp	RA	Resp	Concentration	Name
42:55	4.78e+05	5.35e+05	0.89 y	1.01e+06	116	
44:18	4.53e+05	4.17e+05	1.09 y	8.70e+05	99.8	1,2,3,4,6,7,8-HpCDD

Totals class: Total Tetra-Furans

Entry #: 42

Run: 12 File: 24FEB09M S: 15 I: 1 F: 1  
Acquired: 24-FEB-09 02:31:11

Total Concentration: 351 Unnamed Concentration: 311.338

RT	m1 Resp	m2 Resp RA	Resp	Concentration	Name
23:12	1.88e+04	2.50e+04 0.75 y	4.38e+04	1.44	
23:34	5.38e+04	8.04e+04 0.67 y	1.34e+05	4.41	
23:56	5.67e+04	8.15e+04 0.70 y	1.38e+05	4.54	
24:21	4.05e+05	5.84e+05 0.69 y	9.89e+05	32.5	
24:35	6.26e+05	8.88e+05 0.70 y	1.51e+06	49.7	
24:52	2.07e+05	3.06e+05 0.68 y	5.13e+05	16.9	
25:06	1.57e+05	2.26e+05 0.70 y	3.84e+05	12.6	
25:26	2.57e+05	3.81e+05 0.68 y	6.38e+05	21.0	
25:42	4.61e+05	6.63e+05 0.70 y	1.12e+06	37.0	
26:06	9.39e+05	1.31e+06 0.72 y	2.25e+06	74.0	
26:19	1.18e+05	1.73e+05 0.68 y	2.91e+05	9.56	
26:26	3.69e+04	4.99e+04 0.74 y	8.67e+04	2.85	
26:42	2.83e+05	4.19e+05 0.67 y	7.02e+05	23.1	
26:48	4.93e+05	7.14e+05 0.69 y	1.21e+06	39.7	2,3,7,8-TCDF
26:59	1.56e+04	2.01e+04 0.78 y	3.56e+04	1.17	
27:07	8.80e+04	1.26e+05 0.70 y	2.14e+05	7.03	
27:21	6.85e+04	8.50e+04 0.81 y	1.53e+05	5.04	
28:36	4.22e+04	4.81e+04 0.88 y	9.03e+04	2.97	
28:43	7.09e+04	1.01e+05 0.70 y	1.72e+05	5.66	

Totals class: 1st Fn. Tot Penta-Furans      Entry #: 43

Run: 12      File: 24FEB09M      S: 15 I: 1 F: 1  
Acquired: 24-FEB-09 02:31:11

Total Concentration: 32.6      Unnamed Concentration: 32.600

RT	m1 Resp	m2 Resp	RA	Resp	Concentration	Name
28:35	3.77e+05	2.34e+05	1.61 y	6.10e+05	32.6	

Totals class: Total Penta-Furans

Entry #: 44

Run: 12

File: 24FEB09M

Acquired: 24-FEB-09 02:31:11

S: 15 I: 1 F: 2

Total Concentration: 398

Unnamed Concentration: 259.667

RT	m1 Resp	m2 Resp RA	Resp	Concentration	Name
30:14	1.84e+05	1.09e+05	1.69 y	2.93e+05	
30:22	1.42e+06	9.01e+05	1.57 y	2.32e+06	15.6
31:03	7.01e+04	5.02e+04	1.40 y	1.20e+05	124
31:12	2.17e+05	1.34e+05	1.62 y	3.50e+05	6.42
31:23	2.59e+05	1.79e+05	1.45 y	4.37e+05	18.7
31:39	1.46e+06	9.50e+05	1.53 y	2.41e+06	23.4
32:01	5.43e+05	3.42e+05	1.59 y	8.86e+05	113 1,2,3,7,8-PeCDF
32:57	2.57e+05	1.61e+05	1.59 y	4.18e+05	47.3
33:00	1.86e+05	1.16e+05	1.60 y	3.02e+05	25.7 2,3,4,7,8-PeCDF
33:52	3.24e+04	2.27e+04	1.43 y	5.50e+04	16.1
34:21	6.11e+04	4.05e+04	1.51 y	1.02e+05	2.94
					5.43

Totals class: Total Hexa-Furans

Entry #: 45

Run: 12 File: 24FEB09M

S: 15 I: 1 F: 3

Acquired: 24-FEB-09 02:31:11

Total Concentration: 569

Unnamed Concentration: 159.800

RT	m1 Resp	m2 Resp RA	Resp	Concentration	Name
35:44	7.24e+05	5.67e+05 1.28 y	1.29e+06	83.5	
35:55	2.18e+04	1.54e+04 1.41 y	3.73e+04	2.41	
36:02	2.05e+04	1.68e+04 1.22 y	3.73e+04	2.41	
36:16	4.80e+04	3.79e+04 1.27 y	8.59e+04	5.55	
36:33	2.95e+05	2.35e+05 1.25 y	5.30e+05	34.3	
37:09	9.59e+04	7.74e+04 1.24 y	1.73e+05	11.2	
37:20	2.21e+06	1.74e+06 1.27 y	3.96e+06	247	1,2,3,4,7,8-HxCDF
37:31	9.38e+05	7.40e+05 1.27 y	1.68e+06	110	1,2,3,6,7,8-HxCDF
37:46	6.26e+04	4.85e+04 1.29 y	1.11e+05	7.18	
37:59	3.25e+04	2.90e+04 1.12 y	6.15e+04	3.98	
38:09	7.71e+04	6.71e+04 1.15 y	1.44e+05	9.33	
38:28	1.77e+05	1.45e+05 1.23 y	3.22e+05	17.6	2,3,4,6,7,8-HxCDF
39:57	2.33e+05	1.85e+05 1.26 y	4.18e+05	34.1	1,2,3,7,8,9-HxCDF

Totals class: Total Hepta-Furans

Entry #: 46

Run: 12 File: 24FEB09M S: 15 I: 1 F: 4  
Acquired: 24-FEB-09 02:31:11

Total Concentration: 777 Unnamed Concentration: 170.208

RT	m1 Resp	m2 Resp RA	Resp	Concentration	Name
42:24	2.56e+06	2.49e+06 1.03 y	5.04e+06	475	1,2,3,4,6,7,8-HpCDF
42:55	1.69e+05	1.73e+05 0.98 y	3.42e+05	32.2	
43:13	7.37e+05	7.29e+05 1.01 y	1.47e+06	138	
45:13	7.13e+05	6.83e+05 1.04 y	1.40e+06	133	1,2,3,4,7,8,9-HpCDF

FAL ID: 5340-002-0001-SA      Filename: 04MAR09A      Sam:4      Acquired: 4-MAR-09 11:31:50      iCal: TCDFFAL1-1-13-09  
Client ID: E3YDO      ConCal: ST030409A1 EndCal: ST030409A2  
Results:      GC Column: DB225      Amount: 1.010

Name	Resp	RA	RT	RRF	Conc	Qual	Fac	Noise	DL	#Hom	Rec
2,3,7,8-TCDF	2.29e+06	0.78	y	19:27	1.03	47.7	2.50	-	-	1	98.1
13C-2,3,7,8-TCDF	9.25e+07	0.80	y	19:26	0.87	1940					
13C-1,2,3,4-TCDF	1.08e+08	0.81	y	16:53	-	57.6					

Analyst: 8 Date: 3/4/09

2DF - FORM II-HR CDD  
CDD/CDF TOTAL HOMOLOGUE CONCENTRATION SUMMARY  
HIGH RESOLUTION

SAMPLE No.  
E3YC9

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: SDG NO.: SDG E3YC9

MATRIX: (SOIL/WATER/ASH/TISSUE/OIL) SOIL LAB SAMPLE ID: 5340-001-0001-SA

SAMPLE wt/vol: 1.000 (g/mL): g LAB FILE ID: 24FEB09M Sam: 14

WATER SAMPLE PREP: SPE (SEPF/SPE) DATE RECEIVED: 11-FEB-09

CONCENTRATED EXTRACT VOLUME: 20 (uL) DATE EXTRACTED: 18-FEB-09

INJECTION VOLUME: 2 (uL) % SOLIDS/LIPIDS: 75.15 DATE ANALYZED: 24-FEB-09

GC COLUMN: DB5 ID: 0.25 (mm) DILUTION FACTOR: 1

CONCENTRATION UNITS: (pg/L or ng/kg) ng/kg

HOMOLOGUE DIOXINS	PEAKS	CONCENTRATION	Q	EMPC/EDL
Total TCDD	6	12.6		*
Total PeCDD	2	25.9		*
Total HxCDD	5	*	X	254
Total HpCDD	2	283		*
FURANS				
Total TCDF	18	494		*
Total PeCDF	15	635		*
Total HxCDF	15	822		*
Total HpCDF	4	1060		*

NOTE: Concentrations, Estimated Maximum Possible Concentrations (EMPCs), and Estimated Detection Limits (EDLs) for solid samples are calculated on a dry weight basis (except tissues, which are reported on a wet weight basis with % Lipids). The total homologue concentrations do not affect the TEF (Toxicity Equivalent Factor) calculations.

ANALYST: SC

DATE: 3/2/09

2DF - FORM II-HR CDD  
CDD/CDF TOTAL HOMOLOGUE CONCENTRATION SUMMARY  
HIGH RESOLUTION

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: NA SDG NO.: SDG E3YC9

MATRIX: (SOIL/WATER/ASH/TISSUE/OIL) SOIL LAB SAMPLE ID: 5340-002-0001-SA

SAMPLE wt/vol: 1.010 (g/mL): g LAB FILE ID: 24FEB09M Sam: 15

WATER SAMPLE PREP: SPE (SEPF/SPE) DATE RECEIVED: 11-FEB-09

CONCENTRATED EXTRACT VOLUME: 20 (uL) DATE EXTRACTED: 18-FEB-09

INJECTION VOLUME: 2 (uL) % SOLIDS/LIPIDS: 75.15 DATE ANALYZED: 24-FEB-09

GC COLUMN: DB5 ID: 0.25 (mm) DILUTION FACTOR: 1

CONCENTRATION UNITS: (pg/L or ng/kg) ng/kg

HOMOLOGUE	PEAKS	CONCENTRATION	Q	EMPC/EDL
<b>DIOXINS</b>				
Total TCDD	4	7.61		*
Total PeCDD	1	*	X	12.1
Total HxCDD	5	*	X	209
Total HpCDD	2	216		*
<b>FURANS</b>				
Total TCDF	19	351		*
Total PeCDF	12	431		*
Total HxCDF	13	569		*
Total HpCDF	4	777		*

NOTE: Concentrations, Estimated Maximum Possible Concentrations (EMPCs), and Estimated Detection Limits (EDLs) for solid samples are calculated on a dry weight basis (except tissues, which are reported on a wet weight basis with % Lipids). The total homologue concentrations do not affect the TEF (Toxicity Equivalent Factor) calculations.

ANALYST: SC

DATE: 3/2/09

## CDD/CDF COMPLETE SDG FILE (CSF) INVENTORY SHEET

LABORATORY NAME Frontier Analytical LaboratoryCITY/STATE El Dorado Hills CACASE NO EPADE-2001 SDG NO SDG ESY19 SDG NOS. TO FOLLOW \_\_\_\_\_

TASK ORDER NO. \_\_\_\_\_

CONTRACT NO. EP08W001564

SOW NO. \_\_\_\_\_

All documents delivered in the Complete SDG File must be original documents where possible.  
 (Reference - Exhibit B Section 2.6)


	<u>PAGE NOS.</u>		<u>CHECK</u>	
	<u>FROM</u>	<u>TO</u>	<u>LAB</u>	<u>EPA</u>
1. <u>Inventory Sheet</u> (DC-2) (Do not number)	<u>NA</u>	<u>NA</u>	<u>/</u>	<u> </u>
2. <u>SDG Narrative</u>	<u>1</u>	<u>4</u>	<u>✓</u>	<u> </u>
3. <u>Traffic Report</u>	<u>5</u>	<u>9</u>	<u>✓</u>	<u> </u>
4. <u>CDD/CDF Data</u>				
a. Sample Data				
Sample Data Summary (FORM I-HR CDD-1)	<u>10</u>	<u>11</u>	<u>/</u>	<u> </u>
Toxicity Equivalence Summary (FORM I-HR CDD-2)	<u>12</u>	<u>13</u>	<u>✓</u>	<u> </u>
Second Column confirmation Summary (FORM I-HR CDD-3)	<u>14</u>	<u>15</u>	<u>/</u>	<u> </u>
Selected Ion Current Profile (SICP) for each sample	<u>16</u>	<u>49</u>	<u>/</u>	<u> </u>
Quantitation Reports and Area Summaries	<u>50</u>	<u>71</u>	<u>✓</u>	<u> </u>
Total Homologue Concentration Summary (FORM II-HR CDD)	<u>72</u>	<u>73</u>	<u>✓</u>	<u> </u>
b. Quality Control Data				
Lab Control Sample Summary (FORM III-HR CDD-1)	<u>74</u>	<u>74</u>	<u>✓</u>	<u> </u>
Lab Control Sample Duplicate Summary (FORM III-HR CDD-2)	<u>NA</u>	<u>NA</u>	<u>✓</u>	<u> </u>
Method Blank Summary (FORM IV-HR CDD)	<u>75</u>	<u>75</u>	<u>/</u>	<u> </u>
Window Defining Mix Summary (FORM V-HR CDD-1)	<u>76</u>	<u>79</u>	<u>✓</u>	<u> </u>
Chromatographic Resolution Summary (FORM V-HR CDD-2)	<u>80</u>	<u>83</u>	<u>✓</u>	<u> </u>
Analytical Sequence Summary (FORM V-HR CDD-3)	<u>84</u>	<u>85</u>	<u>/</u>	<u> </u>
c. Calibration Data				
Initial Calibration Data (FORM VI-HR CDD-1 and FORM VI-HR CDD-2), PFK mass resolution, CDD/CDF standard(s) SICPs, Quantitation Reports, and Area Summaries for the initial (five-point) calibration	<u>86</u>	<u>189</u>	<u>/</u>	<u> </u>
Continuing Calibration Data (FORM VII-HR CDD-1 and FORM VII-HR CDD-2), PFK mass resolution, SICPs, Quantitation Reports, and Area Summaries	<u>190</u>	<u>241</u>	<u>✓</u>	<u> </u>
d. Raw Quality Control Data				
Blank Data FORM I-HR CDD-1, CDD-2, CDD-3 (if applicable)	<u>242</u>	<u>243</u>	<u>✓</u>	<u> </u>
Blank Data including SICPs, Quantitation Reports, and Area Summaries for each blank analyzed	<u>244</u>	<u>256-256</u>	<u>✓</u>	<u> </u>
LCS FORM I-HR CDD-1 and CDD-2	<u>258-257</u>	<u>257</u>	<u>/</u>	<u> </u>

	PAGE NOS.	CHECK		
	FROM	TO	LAB	EPA
LCS Data including SICPs, Quantitation Reports, and Area Summaries	<u>258</u>	<u>261</u>	<u>✓</u>	—
<b>5. Miscellaneous Data</b>				
Original preparation and analysis forms or copies of preparation and analysis logbook pages	<u>270</u>	<u>273</u>	<u>✓</u>	—
Internal sample and sample extract transfer Chain of Custody Records	<u>274</u>	<u>274</u>	<u>✓</u>	—
Screening records	<u>NA</u>	<u>NA</u>	<u>✓</u>	—
All instrument output, including strip charts from screening activities (describe or list)				
	<u>NA</u>	<u>NA</u>	<u>✓</u>	—
<b>6. EPA Shipping/Receiving Documents</b>				
Airbills (No. of shipments <u>1</u> )	<u>275</u>	<u>275</u>	<u>✓</u>	—
Chain of Custody Records	<u>276</u>	<u>276</u>	<u>✓</u>	—
Sample Tags	<u>NA</u>	<u>NA</u>	<u>✓</u>	—
Sample Log-In Sheet (Lab & DC-1)	<u>277</u>	<u>277</u>	<u>✓</u>	—
Traffic Report Cover Sheet	<u>NA</u>	<u>NA</u>	<u>✓</u>	—
Miscellaneous Shipping/Receiving Records (describe or list)				
	<u>NA</u>	<u>NA</u>	<u>✓</u>	—
<b>7. Internal Lab Sample Transfer Records and Tracking Sheets</b>				
(Describe or list)				
	<u>NA</u>	<u>NA</u>	<u>✓</u>	—
<b>8. Other Records</b> (describe or list)				
Telephone Communication Log <u>Email:</u>				
	<u>278</u>	<u>278</u>	<u>✓</u>	—

**9. Comments:**

Completed by:

(CLP Lab)



(Signature)

Brian Silversky, Director of Operations

(Print Name & Title)

3/5/05

(Date)

Audited by:

(USEPA)

(Signature)

(Print Name & Title)

(Date)



February 11, 2009

**FAL Project ID: 5340**

Ms. Terese A Van Donsel  
USEPA Region 5  
77 W. Jackson Blvd.  
Chicago IL, 60604

Dear Ms. Donsel,

This letter confirms the receipt of two soil samples at Frontier Analytical Laboratory on February 11, 2009. We received two soil samples and a total of 1 bottle. The bottle arrived in good condition. The samples will be extracted and analyzed using EPA Method DLM02.0 for tetra through octa chlorinated dibenzo dioxins and furans. These samples are associated with EPA Order/Contract No.: **EP08W001564**. This project has a turnaround time of fourteen calendar days at \$900/sample. The following e-mail is a PDF copy of the completed Chain-of-Custody/Traffic Report. The samples have been assigned a Frontier Analytical Laboratory project number **5340** and SDG Number **SDG E3YC9** In order to expedite inquiries, please reference our project number.

These samples represent the first and last samples in SDG Number **SDG E3YC9**.

If you have any questions regarding this project, please feel free to contact us at your convenience.

Thank you for choosing Frontier Analytical Laboratory for your analytical testing needs.

Sincerely,

A handwritten signature in black ink, appearing to read "Benj Miller".

Benjamin Miller



**USEPA Organic Traffic Laboratory Program**  
Part E-PROSW001509  
Analytic Under Part E-PROSW001509

2. Date Shipped: 2/16/2009		3. Chain of Custody Record		Sampler Signature: <i>J. M. Burch</i>		4. For Lab Use Only		
Carrier Name: FedEx	Relinquished By: (Date/Time)	Received By: (Date/Time)				Lab Contract No: EP08W001504	DAS No.: SDG	
Airbill: 8637 05044404	1) <i>MJZ</i> - 2/16/09 16:00	2) <i>Stylusitive Shlifer 045AM</i>				Unit Price: \$700	SDG No.: E3YCA	
Shipped To: Frontier Analytical Lab 5172 Hilldale Circle El Dorado Hills, CA 95762	3)	4)				Transfer To:		
5. ORGANIC SAMPLE NO.	6. MATRIX/ SAMPLER	7. TYPE	8. ANALYSIS/ TURNAROUND	9. TAG No./ PRESERVATIVE/ Bottles	10. STATION LOCATION	11. SAMPLE COLLECT DATE/TIME	12. INORGANIC SAMPLE NO.	13. FOR LAB USE ONLY Sample Condition On Receipt
E3YCA	Soil / Plant extract	Farms	Dioxins / Furans	1/13/2009	January 13, 2009	13:30	5340-001-SA First	good, 2/11/09, 10:45
E3YDP	Soil / Product	GPPS	Dioxins / Furans	1/13/2009	1/13/2009	13:30	5340-002-SA Last	good, 2/11/09, 10:45
14. Shipment for Case Complete?	15. Sample(s) to be used for laboratory QC:	16. Additional Sampler Signature(s):		17. Cooler Temperature Upon Receipt:		18. Chain of Custody Seal Number:		
✓	E3YCA			8°C				
19. Analysis Key:	Type: Comp, Grab (from Box 7)			20. Custody Seal Intact? ✓		21. Shipment Iced? ✓		



## Frontier Analytical Laboratory

### Sample Tracking Log

FAL Project ID: **5340**

Received on: **02/11/2009**

Project Due: **02/26/2009** Storage: **R1**

FAL Sample ID	Dup	Client Project ID	Client Sample ID	Requested Method	Matrix	Sampling Date	Sampling Time
5340-001-SA	0	EP08W001564	E3YC9	DLM02.0 D/F	Soil	01/13/2009	01:30 pm
5340-002-SA	0	EP08W001564	E3YD0	DLM02.0 D/F	Soil	01/13/2009	01:30 pm

000007 of 000279

## Frontier Analytical Laboratory

### Sample Login Form

FAL Project ID: **5340**

Client:	USEPA Region 5
Client Project ID:	EP08W001564
Date Received:	02/11/2009
Time Received:	10:45 am
Received By:	TC
Logged In By:	BM
# of Samples Received:	2
Duplicates:	0
Storage Location:	R1

Method of Delivery:	Fed-Ex
Tracking Number:	868144787343
Shipping Container Received Intact	Yes
Custody seals(s) present?	Yes
Custody seals(s) intact?	Yes
Sample Arrival Temperature (C)	18
Cooling Method	Other
Chain Of Custody Present?	Yes
Return Shipping Container To Client	No
Test for residual Chlorine	No
Thiosulfate Added	No
Earliest Sample Hold Time Expiration	02/18/2009
Adequate Sample Volume	Yes
Anomalies or additional comments:	

Please note that the sample was received at 18 degrees C which is outside the method recommended temperature range of 0 – 4 degrees C. We will proceed with analysis unless instructed otherwise.



USEPA Contract Laboratory Program  
Organic Traffic Report & Chain of Custody Record

Analysis under P# EP08W001564


5340  
18°C

1. Case No.: P07EP08W001564

DAS No.:

SDG No.: SDG E3YC9

2. Date Shipped: 2/16/2009			3. Chain of Custody Record			4. For Lab Use Only		
Carrier Name: FedEx			Relinquished By: (Date/Time)			Received By: (Date/Time)		
Airbill: 8632.05044404			1) <i>Chay</i> 2/16/09 1600			2) <i>Tony Chaitow</i> 2/16/09 10:45AM		
Shipped To: Frontier Analytical Lab 5172 Hillsdale Circle El Dorado Hills, CA 95762			3)			4)		
5. ORGANIC SAMPLE No.	6. MATRIX/ SAMPLER	7. TYPE	8. ANALYSIS/ TURNAROUND	9. TAG No./ PRESERVATIVE/Bottles	10. STATION LOCATION	11. SAMPLE COLLECT DATE/TIME	12. INORGANIC SAMPLE No.	13. FOR LAB USE ONLY Sample Condition On Receipt
E3YC9	SIL / Product	GRAB	Dinxins / Furans	1808807-1m		January 13, 2009 1330		good, 2/11/09, 10:45 5340-001-SA First
E3YD4	SIL / Product	GRAB	Dinxins / Furans	+ 808807-1m Please run a duplicate from same jar		1/13/2009 1330		good, 2/11/09, 10:45 5340-002-SA Last
14. Shipment for Case Complete? <input checked="" type="checkbox"/>			15. Sample(s) to be used for laboratory analysis: E3YC9			16. Cooler Temperature on Receipt: 18°C		
17. Analysis Key: Type: Comp, Grab (from Box 7)			18. Chain of Custody Seal Number: 5340-002-SA			19. Custody Seal Intact? <input checked="" type="checkbox"/> 20. Shipment Iced? <input type="checkbox"/>		
<p>22. TR Number: PR provides preliminary results. Requests for preliminary results will be sent to: Sample Management Office, Attn: Heather Bauer, 1501</p> <p>23. Laboratory COPY Page 1 of 6</p>								



3DFA - FORM III-HR CDD  
CDD/CDF LAB CONTROL SAMPLE SUMMARY  
HIGH RESOLUTION

SAMPLE No.  
DLCS01

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: SDG NO.: SDG E3YC9

MATRIX: (SOIL/WATER/ASH/TISSUE/OIL) SOIL LAB SAMPLE ID: 1652-001-0001-OPR

SAMPLE wt/vol: 1.000 (g/mL): g LAB FILE ID: 24FEB09M Sam: 12

WATER SAMPLE PREP: SPE (SEPF/SPE) DATE RECEIVED: 18-FEB-09

CONCENTRATED EXTRACT VOLUME: 20 (uL) DATE EXTRACTED: 18-FEB-09

INJECTION VOLUME: 2 (uL) DATE ANALYZED: 23-FEB-09

GC COLUMN: DB5 ID: 0.25 (mm) DILUTION FACTOR: 1

CONCENTRATION UNITS: (pg/L or ng/kg) ng/kg

SPIKE ANALYTE	SPIKE ADDED	AMOUNT RECOVERED	PERCENT RECOVERY	QC #	QC LIMITS
2,3,7,8-Tcdd	10	9.67	96.7		67-158
2,3,7,8-Tcdf	10	8.38	83.8		75-158
1,2,3,7,8-PeCdf	50	46.1	92.3		80-134
1,2,3,7,8-PeCdd	50	50.1	100		70-142
2,3,4,7,8-PeCdf	50	46.4	92.7		68-160
1,2,3,4,7,8-HxCdf	50	51.1	102		72-134
1,2,3,6,7,8-HxCdf	50	50.5	101		84-130
1,2,3,4,7,8-HxCdd	50	47.0	94.1		70-164
1,2,3,6,7,8-HxCdd	50	46.0	92.1		76-134
1,2,3,7,8,9-HxCdd	50	46.2	92.4		64-162
2,3,4,6,7,8-HxCdf	50	50.1	100		70-156
1,2,3,7,8,9-HxCdf	50	50.4	101		78-130
1,2,3,4,6,7,8-HpCDF	50	50.8	102		82-132
1,2,3,4,6,7,8-HpCdd	50	43.4	86.8		70-140
1,2,3,4,7,8,9-HpCDF	50	51.8	104		78-138
OCDD	100	89.2	89.2		78-144
OCDF	100	94.3	94.3		63-170

# Column to be used to flag values outside Quality Control (QC) Limits.

Laboratory Control Sample Recovery: 0 Outside limits of 17 total.

ANALYST: R

DATE: 3/2/09

4DF - FORM IV-HR CDD  
CDD/CDF METHOD BLANK SUMMARY  
HIGH RESOLUTION

SAMPLE No.  
DFBLK01

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: SDG NO.: SDG E3YC9

MATRIX: (SOIL/WATER/ASH/TISSUE/OIL) SOIL LAB SAMPLE ID: 1652-001-0001-MB

WATER SAMPLE PREP: SPE (SEPF/SPE) LAB FILE ID: 24FEB09M Sam: 13

GC COLUMN: DB5 ID: 0.25 (mm) DATE EXTRACTED: 18-FEB-09

INSTRUMENT ID: FAL3 DATE ANALYZED: 24-FEB-09

THIS METHOD BLANK APPLIES TO LABORATORY CONTROL SAMPLES (LCSS).

EPA SAMPLE NO.	LAB SAMPLE ID	LAB FILE ID	DATE ANALYZED
OPR	1652-001-0001-OPR	24FEB09M 12	23-FEB-09
Method Blank	1652-001-0001-MB	24FEB09M 13	24-FEB-09
E3YC9	5340-001-0001-SA	24FEB09M 14	24-FEB-09
E3YD0	5340-002-0001-SA	24FEB09M 15	24-FEB-09

ANALYST: 

DATE:  3/2/09

SDFA - FORM V-HR CDD-1  
CDD/CDF WINDOW DEFINING MIX (WDM) SUMMARY  
HIGH RESOLUTION

SAMPLE No.  
1613 CS3 (080827J)

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: SDG NO.: SDG E3YC9

GC COLUMN: DB5 ID: 0.25 (mm) LAB FILE ID: 24FEB09M Sam: 11

INSTRUMENT ID: FAL3 DATE ANALYZED: 23-FEB-09

TIME ANALYZED: 22:50:13

CDD/CDF	RT FIRST ELUTING	RT LAST ELUTING
TCDD	24:31	28:28
TCDF	23:10	28:42
PeCDD	30:21	33:55
PeCDF	28:32	34:21
HxCDD	36:14	39:18
HxCDF	35:21	39:52
HpCDD	42:55	44:17
HpCDF	42:24	45:13

ANALYST: [Signature]

DATE: 3/3/09

5DFA - FORM V-HR CDD-1  
CDD/CDF WINDOW DEFINING MIX (WDM) SUMMARY  
HIGH RESOLUTION

SAMPLE No.  
1613 CS3 (080827J)

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: SDG NO.: SDG E3YC9

GC COLUMN: DB5 ID: 0.25 (mm) LAB FILE ID: 24FEB09M Sam: 18

INSTRUMENT ID: FAL3 DATE ANALYZED: 24-FEB-09

TIME ANALYZED: 05:17:02

CDD/CDF	RT FIRST ELUTING	RT LAST ELUTING
TCDD	24:30	28:27
TCDF	23:09	28:41
PeCDD	30:21	33:55
PeCDF	28:31	34:19
HxCDD	36:13	39:18
HxCDF	35:20	39:53
HpCDD	42:54	44:17
HpCDF	42:22	45:13

ANALYST: T

DATE: 3/2/09

SDFA - FORM V-HR CDD-1  
TCDF WINDOW DEFINING MIX (WDM) SUMMARY  
HIGH RESOLUTION

SAMPLE No.  
1613 CS3 (080827J)

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: SDG NO.: SDG E3YC9

GC COLUMN: DB225 ID: 0.25 (mm) LAB FILE ID: 04MAR09A Sam: 1

INSTRUMENT ID: FAL1 DATE ANALYZED: 4-MAR-09

TIME ANALYZED: 09:46:39

CDD/CDF	RT FIRST	RT LAST
ELUTING	ELUTING	
TCDF	13:37	22:07

ANALYST: J DATE: 3/4/09

SDFA - FORM V-HR CDD-1  
TCDF WINDOW DEFINING MIX (WDM) SUMMARY  
HIGH RESOLUTION

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: SDG NO.: SDG E3YC9

GC COLUMN: DB225 ID: 0.25 (mm) LAB FILE ID: 04MAR09A Sam: 5

INSTRUMENT ID: FAL1 DATE ANALYZED: 4-MAR-09

TIME ANALYZED: 12:09:12

CDD/CDF	RT FIRST	RT LAST
	ELUTING	ELUTING

TCDF	13:31	22:02
------	-------	-------

ANALYST: 

DATE: 3/4/09

5DFB - FORM V-HR CDD-2  
CDD/CDF CHROMATOGRAPHIC RESOLUTION SUMMARY  
HIGH RESOLUTION

SAMPLE No.  
1613 CS3 (080827J)

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: SDG NO.: SDG E3YC9

GC COLUMN: DB5 ID: 0.25 (mm) LAB FILE ID: 24FEB09M Sam: 11

INSTRUMENT ID: FAL3 DATE ANALYZED: 23-FEB-09

TIME ANALYZED: 22:50:13

Percent Valley determination for DB-5 (or equivalent) column-  
For the column performance solution beginning the 12-hour period:

1238-TCDD/2378-TCDD: 15.2 %

QUALITY CONTROL (QC) LIMITS:

Percent Valley between the TCDD isomers must be less than or equal to 25%.

Percent Valley determination for DB-225 (or equivalent) column-  
For the column performance solution beginning the 12-hour period:

2347-TCDF/2378-TCDF: NA

QC Limits:

Percent Valley between the TCDD/TCDF isomers must be less than or equal to 25%.

ANALYST: SC

DATE: 3/2/09

5DFB - FORM V-HR CDD-2  
CDD/CDF CHROMATOGRAPHIC RESOLUTION SUMMARY  
HIGH RESOLUTION

SAMPLE No.  
1613 CS3 (080827J)

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: SDG NO.: SDG E3YC9

GC COLUMN: DB5 ID: 0.25 (mm) LAB FILE ID: 24FEB09M Sam: 18

INSTRUMENT ID: FAL3 DATE ANALYZED: 24-FEB-09

TIME ANALYZED: 05:17:02

Percent Valley determination for DB-5 (or equivalent) column-  
For the column performance solution beginning the 12-hour period:

1238-TCDD/2378-TCDD: 13.8 %

QUALITY CONTROL (QC) LIMITS:

Percent Valley between the TCDD isomers must be less than or equal to 25%.

Percent Valley determination for DB-225 (or equivalent) column-  
For the column performance solution beginning the 12-hour period:

2347-TCDF/2378-TCDF: NA

QC Limits:

Percent Valley between the TCDD/TCDF isomers must be less than or equal to 25%.

ANALYST: E

DATE: 3/2/09

5DFB - FORM V-HR CDD-2  
TCDF CHROMATOGRAPHIC RESOLUTION SUMMARY  
HIGH RESOLUTION

SAMPLE No.  
1613 CS3 (080827J)

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: SDG NO.: SDG E3YC9

GC COLUMN: DB225 ID: 0.25 (mm) LAB FILE ID: 04MAR09A Sam: 1

INSTRUMENT ID: FAL1 DATE ANALYZED: 4-MAR-09

TIME ANALYZED: 09:46:39

Percent Valley determination for DB-5 (or equivalent) column-  
For the column performance solution beginning the 12-hour period:

1238-TCDD/2378-TCDD: NA

QUALITY CONTROL (QC) LIMITS:

Percent Valley between the TCDD isomers must be less than or equal to 25%.

Percent Valley determination for DB-225 (or equivalent) column-  
For the column performance solution beginning the 12-hour period:

2347-TCDF/2378-TCDF: 7.25 %

QC Limits:

Percent Valley between the TCDD/TCDF isomers must be less than or equal to 25%.

ANALYST: J DATE: 3/4/09

5DFB - FORM V-HR CDD-2  
TCDF CHROMATOGRAPHIC RESOLUTION SUMMARY  
HIGH RESOLUTION

SAMPLE No.  
1613 CS3 (080827J)

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: SDG NO.: SDG E3YC9

GC COLUMN: DB225 ID: 0.25 (mm) LAB FILE ID: 04MAR09A Sam: 5

INSTRUMENT ID: FAL1 DATE ANALYZED: 4-MAR-09

TIME ANALYZED: 12:09:12

Percent Valley determination for DB-5 (or equivalent) column-  
For the column performance solution beginning the 12-hour period:

1238-TCDD/2378-TCDD: NA

QUALITY CONTROL (QC) LIMITS:

Percent Valley between the TCDD isomers must be less than or equal to 25%.

Percent Valley determination for DB-225 (or equivalent) column-  
For the column performance solution beginning the 12-hour period:

2347-TCDF/2378-TCDF: 5.80 %

QC Limits:

Percent Valley between the TCDD/TCDF isomers must be less than or equal to 25%.

ANALYST: J

DATE: 3/4/09

5DFC - FORM V-HR CDD-3  
CDD/CDF ANALYTICAL SEQUENCE SUMMARY  
HIGH RESOLUTION

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: SDG NO.: SDG E3YC9

GC COLUMN: DB5 ID: 0.25 (mm) INSTRUMENT ID: FAL3

INIT. CALIB. DATE (S) : 5-FEB-09

INIT. CALIB. TIMES: 16:15:11, 17:10:30, 18:05:53, 15:20:00, 19:01:11, 19:56:22

THE ANALYTICAL SEQUENCE OF STANDARDS, SAMPLES, BLANKS, AND LABORATORY CONTROL SAMPLES (LCSS) IS AS FOLLOWS:

EPA SAMPLE NO.	LAB SAMPLE ID	LAB FILE ID	DATE ANALYZED	TIME ANALYZED
ZZZZ	ZZZZ	24FEB09M	23-FEB-09	13:37:17
ZZZZ	ZZZZ	24FEB09M	23-FEB-09	14:32:35
ZZZZ	ZZZZ	24FEB09M	23-FEB-09	15:27:54
ZZZZ	ZZZZ	24FEB09M	23-FEB-09	16:23:13
ZZZZ	ZZZZ	24FEB09M	23-FEB-09	17:18:27
ZZZZ	ZZZZ	24FEB09M	23-FEB-09	18:13:42
ZZZZ	ZZZZ	24FEB09M	23-FEB-09	19:09:01
ZZZZ	ZZZZ	24FEB09M	23-FEB-09	20:04:20
ZZZZ	ZZZZ	24FEB09M	23-FEB-09	20:59:39
ZZZZ	ZZZZ	24FEB09M	23-FEB-09	21:54:57
1613 CS3 (080827J)	ST024109M3	24FEB09M	23-FEB-09	22:50:13
OPR	1652-001-0001-OPR	24FEB09M	23-FEB-09	23:45:23
Method Blank	1652-001-0001-MB	24FEB09M	24-FEB-09	00:40:37
E3YC9	5340-001-0001-SA	24FEB09M	24-FEB-09	01:35:56
E3YD0	5340-002-0001-SA	24FEB09M	24-FEB-09	02:31:11
ZZZZ	ZZZZ	24FEB09M	24-FEB-09	03:26:28
ZZZZ	ZZZZ	24FEB09M	24-FEB-09	04:21:47
1613 CS3 (080827J)	ST022409M4	24FEB09M	24-FEB-09	05:17:02

ANALYST: S

DATE: 3/2/09

5DFC - FORM V-HR CDD-3  
TCDF ANALYTICAL SEQUENCE SUMMARY  
HIGH RESOLUTION

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: SDG NO.: SDG E3YC9

GC Column: DB225 ID: 0.25 (mm) INSTRUMENT ID: FAL1

INIT. CALIB. DATE (S) : 13-JAN-09

INIT. CALIB. TIMES: 11:41:31, 12:16:35, 11:06:28, 12:51:39, 13:26:43

THE ANALYTICAL SEQUENCE OF STANDARDS, SAMPLES, BLANKS, AND LABORATORY CONTROL SAMPLES (LCSS) IS AS FOLLOWS:

EPA SAMPLE NO.	LAB SAMPLE ID	LAB FILE ID	DATE ANALYZED	TIME ANALYZED
1613 CS3 (080827J)	ST030409A1	04MAR09A	4-MAR-09	09:46:39
ZZZZ	ZZZZ	04MAR09A	4-MAR-09	10:21:42
E3YC9	5340-001-0001-SA	04MAR09A	4-MAR-09	10:56:46
E3YD0	5340-002-0001-SA	04MAR09A	4-MAR-09	11:31:50
1613 CS3 (080827J)	ST030409A2	04MAR09A	4-MAR-09	12:09:12

ANALYST: J DATE: 3/4/09

6DFA - FORM VI-HR CDD-1  
CDD/CDF INITIAL CALIBRATION RESPONSE FACTOR SUMMARY  
HIGH RESOLUTION

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: NA SDG NO.: SDG E3YC9

GC COLUMN: DB5 ID: 0.25 (mm) INSTRUMENT ID: FAL3

INIT. CALIB. DATE(S): 5-FEB-09 5-FEB-09 5-FEB-09 5-FEB-09 5-FEB-09 5-FEB-09


INIT. CALIB. TIMES: 16:15:11 17:10:30 18:05:53 15:20:00 19:01:11 19:56:22

TARGET ANALYTES	RR/RRF						MEAN RR/RRF	%RSD	QC LIMITS
	CS0	CS1	CS2	CS3	CS4	CS5			
2,3,7,8-TCDD	1.06	1.08	0.98	0.98	1.04	1.10	1.04	4.83	+/-20%
2,3,7,8-TCDF	1.51	1.43	1.29	1.24	1.25	1.31	1.34	8.09	+/-20%
1,2,3,7,8-PeCDF	1.01	0.97	1.01	1.03	1.08	1.08	1.03	4.26	+/-20%
1,2,3,7,8-PeCDD	0.90	0.90	0.87	0.88	0.92	0.94	0.90	3.00	+/-20%
2,3,4,7,8-PeCDF	0.84	0.76	0.80	0.77	0.83	0.86	0.81	4.85	+/-20%
1,2,3,4,7,8-HxCDF	1.29	1.29	1.30	1.33	1.44	1.47	1.35	5.94	+/-20%
1,2,3,6,7,8-HxCDF	1.18	1.22	1.21	1.23	1.31	1.35	1.25	5.32	+/-20%
1,2,3,4,7,8-HxCDD	1.45	1.42	1.40	1.38	1.43	1.49	1.43	2.75	+/-20%
1,2,3,6,7,8-HxCDD	1.04	1.04	1.00	0.97	1.03	1.04	1.02	2.92	+/-20%
1,2,3,7,8,9-HxCDD (1)	1.13	1.07	1.02	1.02	1.01	1.06	1.05	4.39	+/-20%
2,3,4,6,7,8-HxCDF	1.46	1.50	1.46	1.51	1.60	1.65	1.53	5.12	+/-20%
1,2,3,7,8,9-HxCDF	1.17	1.15	1.13	1.14	1.24	1.28	1.19	5.21	+/-20%
1,2,3,4,6,7,8-HpCDF	1.04	1.07	1.05	1.07	1.16	1.21	1.10	6.18	+/-20%
1,2,3,4,6,7,8-HpCDD	1.03	1.01	0.96	0.93	0.98	1.01	0.99	3.97	+/-20%
1,2,3,4,7,8,9-HpCDF	0.98	1.00	0.97	1.00	1.08	1.12	1.02	6.07	+/-20%
OCDD	0.97	0.94	0.90	0.90	0.95	0.97	0.94	3.39	+/-20%
OCDF	0.87	0.71	0.76	0.70	0.76	0.79	0.76	7.87	+/-20%

LABELED COMPOUNDS

13C-2,3,7,8-TCDD	0.91	0.91	0.94	0.91	0.96	0.98	0.94	3.42	+/-35%
13C-1,2,3,7,8-PeCDD	0.70	0.72	0.71	0.75	0.80	0.86	0.75	8.16	+/-35%
13C-1,2,3,4,7,8-HxCDD	0.99	1.03	1.03	1.03	1.05	1.03	1.03	1.78	+/-35%
13C-1,2,3,6,7,8-HxCDD	1.05	1.06	1.11	1.07	1.09	1.10	1.08	2.30	+/-35%
13C-1,2,3,4,6,7,8-HpCDD	0.90	0.87	0.89	0.93	0.88	0.90	0.89	2.39	+/-35%
13C-OCDD	0.64	0.61	0.61	0.67	0.61	0.67	0.64	4.42	+/-35%
13C-2,3,7,8-TCDF	0.85	0.87	0.85	0.86	0.86	0.88	0.86	1.26	+/-35%
13C-1,2,3,7,8-PeCDF	0.69	0.76	0.71	0.75	0.77	0.87	0.76	8.30	+/-35%
13C-2,3,4,7,8-PeCDF	0.73	0.77	0.72	0.77	0.79	0.87	0.78	6.96	+/-35%
13C-1,2,3,4,7,8-HxCDF	1.35	1.35	1.37	1.36	1.34	1.29	1.34	2.16	+/-35%
13C-1,2,3,6,7,8-HxCDF	1.37	1.39	1.43	1.41	1.43	1.37	1.40	1.93	+/-35%
13C-1,2,3,7,8,9-HxCDF	1.01	0.99	1.02	1.00	0.99	1.02	1.01	1.19	+/-35%
13C-2,3,4,6,7,8-HxCDF	1.26	1.28	1.33	1.28	1.30	1.28	1.29	1.75	+/-35%
13C-1,2,3,4,6,7,8-HpCDF	1.16	1.15	1.18	1.16	1.16	1.13	1.16	1.33	+/-35%
13C-1,2,3,4,7,8,9-HpCDF	0.96	0.94	0.96	0.98	0.93	0.98	0.96	2.06	+/-35%
13C-OCDF	0.93	0.91	0.93	0.97	0.93	1.02	0.95	4.34	+/-35%
37Cl-2,3,7,8-TCDD	1.01	0.80	0.83	0.79	0.88	0.96	0.88	9.92	+/-35%

(1) The Relative Response (RR) is calculated based on the labeled analogs of the other two HxCDDs.

Analyst: 

Date:  3/2/09

6DFB - FORM VI-HR CDD-2  
 CDD/CDF INITIAL CALIBRATION ION ABUNDANCE RATIO SUMMARY  
 HIGH RESOLUTION

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: NA SDG NO.: SDG E3YC9

GC COLUMN: DB5 ID: 0.25 (mm) INSTRUMENT ID: FAL3

INIT. CALIB. DATE(S): 5-FEB-09 5-FEB-09 5-FEB-09 5-FEB-09 5-FEB-09 5-FEB-09

INIT. CALIB. TIMES: 16:15:11 17:10:30 18:05:53 15:20:00 19:01:11 19:56:22

TARGET ANALYTES	SELECTED IONS	ION ABUNDANCE RATIO					FLAG	ION RATIO QC LIMITS
		CS0	CS1	CS2	CS3	CS4		
2,3,7,8-TCDD	320/322	0.75	0.81	0.71	0.76	0.75	0.76	0.65-0.89
2,3,7,8-TCDF	304/306	0.68	0.68	0.71	0.72	0.73	0.73	0.65-0.89
1,2,3,7,8-PeCDF	340/342	1.42	1.44	1.47	1.53	1.51	1.52	1.32-1.78
1,2,3,7,8-PeCDD	356/358	1.77	1.64	1.66	1.63	1.63	1.61	1.32-1.78
2,3,4,7,8-PeCDF	340/342	1.36	1.38	1.45	1.51	1.53	1.51	1.32-1.78
1,2,3,4,7,8-HxCDF	374/376	1.26	1.25	1.26	1.26	1.25	1.26	1.05-1.43
1,2,3,6,7,8-HxCDF	374/376	1.29	1.28	1.22	1.25	1.25	1.24	1.05-1.43
1,2,3,4,7,8-HxCDD	390/392	1.26	1.36	1.30	1.32	1.30	1.30	1.05-1.43
1,2,3,6,7,8-HxCDD	390/392	1.25	1.25	1.32	1.32	1.31	1.30	1.05-1.43
1,2,3,7,8,9-HxCDD	390/392	1.20	1.33	1.29	1.30	1.32	1.30	1.05-1.43
2,3,4,6,7,8-HxCDF	374/376	1.23	1.24	1.26	1.26	1.25	1.24	1.05-1.43
1,2,3,7,8,9-HxCDF	374/376	1.24	1.28	1.27	1.26	1.25	1.26	1.05-1.43
1,2,3,4,6,7,8-HpCDF	408/410	0.97	1.04	1.03	1.03	1.03	1.03	0.88-1.20
1,2,3,4,6,7,8-HpCDD	424/426	1.08	1.10	1.04	1.03	1.04	1.04	0.88-1.20
1,2,3,4,7,8,9-HpCDF	408/410	0.96	1.00	1.02	1.03	1.05	1.04	0.88-1.20
OCDD	458/460	0.94	0.91	0.91	0.91	0.90	0.90	0.76-1.02
OCDF	442/444	0.85	0.91	0.91	0.90	0.91	0.91	0.76-1.02
<b>LABELED COMPOUNDS</b>								
13C-2,3,7,8-TCDD	332/334	0.79	0.79	0.78	0.78	0.79	0.79	0.65-0.89
13C-1,2,3,7,8-PeCDD	368/370	1.76	1.76	1.78	1.77	1.77	1.76	1.32-1.78
13C-1,2,3,4,7,8-HxCDD	402/404	1.25	1.26	1.28	1.27	1.30	1.26	1.05-1.43
13C-1,2,3,6,7,8-HxCDD	402/404	1.24	1.26	1.25	1.26	1.24	1.25	1.05-1.43
13C-1,2,3,4,6,7,8-HpCDD	436/438	1.11	1.11	1.12	1.10	1.10	1.10	0.88-1.20
13C-OCDD	470/472	0.91	0.90	0.90	0.92	0.91	0.91	0.76-1.02
13C-2,3,7,8-TCDF	316/318	0.82	0.82	0.82	0.81	0.82	0.82	0.65-0.89
13C-1,2,3,7,8-PeCDF	352/354	1.68	1.68	1.70	1.69	1.71	1.68	1.32-1.78
13C-2,3,4,7,8-PeCDF	352/354	1.69	1.67	1.69	1.70	1.71	1.68	1.32-1.78
13C-1,2,3,4,7,8-HxCDF	384/386	0.53	0.53	0.53	0.53	0.53	0.53	0.43-0.59
13C-1,2,3,6,7,8-HxCDF	384/386	0.53	0.53	0.53	0.53	0.54	0.54	0.43-0.59
13C-1,2,3,7,8,9-HxCDF	384/386	0.53	0.53	0.54	0.54	0.53	0.54	0.43-0.59
13C-2,3,4,6,7,8-HxCDF	384/386	0.53	0.54	0.53	0.53	0.53	0.53	0.43-0.59
13C-1,2,3,4,6,7,8-HpCDF	418/420	0.43	0.43	0.42	0.43	0.43	0.43	0.37-0.51
13C-1,2,3,4,7,8,9-HpCDF	418/420	0.42	0.43	0.42	0.42	0.43	0.43	0.37-0.51
13C-OCDF	454/456	0.91	0.91	0.91	0.91	0.92	0.92	0.76-1.02
<b>INTERNAL STANDARDS</b>								
13C-1,2,3,4-TCDD	332/334	0.80	0.79	0.79	0.79	0.79	0.80	0.65-0.89
13C-1,2,3,4-TCDF	316/318	0.82	0.82	0.82	0.82	0.82	0.83	0.65-0.89
13C-1,2,3,7,8,9-HxCDD	402/404	1.25	1.25	1.26	1.27	1.25	1.24	1.05-1.43

Quality Control (QC) limits represent +/-15% window around the theoretical ion abundance ratio. The laboratory must flag any analyte in any calibration solution which does not meet the ion abundance ratio QC limit by placing an asterik in the flag column

Analyst: R

Date: 3/3/09

5DFC - FORM V-HR CDD-3  
CDD/CDF ANALYTICAL SEQUENCE SUMMARY  
HIGH RESOLUTION

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: SDG NO.: SDG E3YC9

GC COLUMN: DB5 ID: 0.25 (mm) INSTRUMENT ID: FAL3

INIT. CALIB. DATE (S) : 5-FEB-09

INIT. CALIB. TIMES: 16:15:11, 17:10:30, 18:05:53, 15:20:00, 19:01:11, 19:56:22

THE ANALYTICAL SEQUENCE OF STANDARDS, SAMPLES, BLANKS, AND LABORATORY CONTROL SAMPLES (LCSS) IS AS FOLLOWS:

EPA SAMPLE NO.	LAB SAMPLE ID	LAB FILE ID	DATE ANALYZED	TIME ANALYZED
1613 CS3 (080827J)	ST020509M3	05FEB09M	5-FEB-09	15:20:00
1613 CS0 (080827J)	ST020509M0	05FEB09M	5-FEB-09	16:15:11
1613 CS1 (080827J)	ST020509M1	05FEB09M	5-FEB-09	17:10:30
1613 CS2 (080827J)	ST020509M2	05FEB09M	5-FEB-09	18:05:53
1613 CS4 (080827J)	ST020509M4	05FEB09M	5-FEB-09	19:01:11
1613 CS5 (080827J)	ST020509M5	05FEB09M	5-FEB-09	19:56:22
ZZZZ	ZZZZ	05FEB09M	5-FEB-09	20:51:38
ZZZZ	ZZZZ	05FEB09M	5-FEB-09	21:47:01
ZZZZ	ZZZZ	05FEB09M	5-FEB-09	22:42:20
ZZZZ	ZZZZ	05FEB09M	5-FEB-09	23:37:31
ZZZZ	ZZZZ	05FEB09M	6-FEB-09	00:32:42
ZZZZ	ZZZZ	05FEB09M	6-FEB-09	01:27:56
ZZZZ	ZZZZ	05FEB09M	6-FEB-09	02:23:11
ZZZZ	ZZZZ	05FEB09M	6-FEB-09	03:18:26
ZZZZ	ZZZZ	05FEB09M	6-FEB-09	04:13:45
ZZZZ	ZZZZ	05FEB09M	6-FEB-09	05:09:07
ZZZZ	ZZZZ	05FEB09M	6-FEB-09	06:04:26
ZZZZ	ZZZZ	05FEB09M	6-FEB-09	06:59:45
ZZZZ	ZZZZ	05FEB09M	6-FEB-09	07:55:04
ZZZZ	ZZZZ	05FEB09M	6-FEB-09	08:50:22
ZZZZ	ZZZZ	05FEB09M	6-FEB-09	09:45:41

ANALYST: SL

DATE: 3/2/09

## Frontier Analytical Laboratory

Data Filename: 05FEB09M

Analyte:	Cal: PCDDFAL3-2-5-09						S1			S5			S6		
	Name	RRF	S. D.	%RSD	RRF#1	RRF#2	RRF#3	RRF#4	RRF#5	RRF#6					
2,3,7,8-TCDD		1.04	0.0502	4.83 %	1.06	1.08	0.98	0.98	1.04	1.10					
1,2,3,7,8-PeCDD		0.90	0.0270	3.00 %	0.90	0.90	0.87	0.88	0.92	0.94					
1,2,3,4,7,8-HxCDD		1.43	0.0393	2.75 %	1.45	1.42	1.40	1.38	1.43	1.49					
1,2,3,6,7,8-HxCDD		1.02	0.0297	2.92 %	1.04	1.04	1.00	0.97	1.03	1.04					
1,2,3,7,8,9-HxCDD		1.05	0.0462	4.39 %	1.13	1.07	1.02	1.02	1.01	1.06					
1,2,3,4,6,7,8-HpCDD		0.99	0.0392	3.97 %	1.03	1.01	0.96	0.93	0.98	1.01					
OCDD		0.94	0.0319	3.39 %	0.97	0.94	0.90	0.90	0.95	0.97					
2,3,7,8-TCDF		1.34	0.108	8.09 %	1.51	1.43	1.29	1.24	1.25	1.31					
1,2,3,7,8-PeCDF		1.03	0.0439	4.26 %	1.01	0.97	1.01	1.03	1.08	1.08					
2,3,4,7,8-PeCDF		0.81	0.0392	4.85 %	0.84	0.76	0.80	0.77	0.83	0.86					
1,2,3,4,7,8-HxCDF		1.35	0.0805	5.94 %	1.29	1.29	1.30	1.33	1.44	1.47					
1,2,3,6,7,8-HxCDF		1.25	0.0665	5.32 %	1.18	1.22	1.21	1.23	1.31	1.35					
2,3,4,6,7,8-HxCDF		1.53	0.0784	5.12 %	1.46	1.50	1.46	1.51	1.60	1.65					
1,2,3,7,8,9-HxCDF		1.19	0.0618	5.21 %	1.17	1.15	1.13	1.14	1.24	1.28					
1,2,3,4,6,7,8-HpCDF		1.10	0.0678	6.18 %	1.04	1.07	1.05	1.07	1.16	1.21					
1,2,3,4,7,8,9-HpCDF		1.02	0.0621	6.07 %	0.98	1.00	0.97	1.00	1.08	1.12					
OCDF		0.76	0.0600	7.87 %	0.87	0.71	0.76	0.70	0.76	0.79					
13C-2,3,7,8-TCDD		0.94	0.0320	3.42 %	0.91	0.91	0.94	0.91	0.96	0.98					
13C-1,2,3,7,8-PeCDD		0.75	0.0615	8.16 %	0.70	0.72	0.71	0.75	0.80	0.86					
13C-1,2,3,4,7,8-HxCDD		1.03	0.0183	1.78 %	0.99	1.03	1.03	1.03	1.05	1.03					
13C-1,2,3,6,7,8-HxCDD		1.08	0.0248	2.30 %	1.05	1.06	1.11	1.07	1.09	1.10					
13C-1,2,3,4,6,7,8-HpCDD		0.89	0.0213	2.39 %	0.90	0.87	0.89	0.93	0.88	0.90					
13C-OCDD		0.64	0.0281	4.42 %	0.64	0.61	0.61	0.67	0.61	0.67					
13C-2,3,7,8-TCDF		0.86	0.0109	1.26 %	0.85	0.87	0.85	0.86	0.86	0.88					
13C-1,2,3,7,8-PeCDF		0.76	0.0629	8.30 %	0.69	0.76	0.71	0.75	0.77	0.87					
13C-2,3,4,7,8-PeCDF		0.78	0.0540	6.96 %	0.73	0.77	0.72	0.77	0.79	0.87					
13C-1,2,3,4,7,8-HxCDF		1.34	0.0290	2.16 %	1.35	1.35	1.37	1.36	1.34	1.29					
13C-1,2,3,6,7,8-HxCDF		1.40	0.0270	1.93 %	1.37	1.39	1.43	1.41	1.43	1.37					
13C-2,3,4,6,7,8-HxCDF		1.29	0.0225	1.75 %	1.26	1.28	1.33	1.28	1.30	1.28					
13C-1,2,3,7,8,9-HxCDF		1.01	0.0120	1.19 %	1.01	0.99	1.02	1.00	0.99	1.02					
13C-1,2,3,4,6,7,8-HpCDF		1.16	0.0154	1.33 %	1.16	1.15	1.18	1.16	1.16	1.13					
13C-1,2,3,4,7,8,9-HpCDF		0.96	0.0197	2.06 %	0.96	0.94	0.96	0.98	0.93	0.98					
13C-OCDF		0.95	0.0411	4.34 %	0.93	0.91	0.93	0.97	0.93	1.02					
37Cl-2,3,7,8-TCDD		0.88	0.0873	9.92 %	1.01	0.80	0.83	0.79	0.88	0.96					
13C-1,2,3,4-TCDD		-	-	- %	-	-	-	-	-	-					
13C-1,2,3,4-TCDF		-	-	- %	-	-	-	-	-	-					
13C-1,2,3,7,8,9-HxCDD		-	-	- %	-	-	-	-	-	-					
Total Tetra-Dioxins		1.04	0.0502	4.83 %	1.06	1.08	0.98	0.98	1.04	1.10					
Total Penta-Dioxins		0.90	0.0270	3.00 %	0.90	0.90	0.87	0.88	0.92	0.94					
Total Hexa-Dioxins		1.16	0.0326	2.81 %	1.20	1.17	1.13	1.12	1.15	1.19					
Total Hepta-Dioxins		0.99	0.0392	3.97 %	1.03	1.01	0.96	0.93	0.98	1.01					
Total Tetra-Furans		1.34	0.108	8.09 %	1.51	1.43	1.29	1.24	1.25	1.31					
1st Fn. Tot Penta-Furans		0.92	0.0382	4.16 %	0.92	0.86	0.90	0.90	0.95	0.97					
Total Penta-Furans		0.92	0.0382	4.16 %	0.92	0.86	0.90	0.90	0.95	0.97					
Total Hexa-Furans		1.34	0.0713	5.34 %	1.28	1.29	1.29	1.31	1.41	1.45					
Total Hepta-Furans		1.06	0.0648	6.09 %	1.01	1.04	1.01	1.04	1.12	1.17					

Analyst: EDate: 2/6/09

Run #1      Filename 05FEB09M      S: 2      Acquired: 5-FEB-09 16:15:11      Cal: PCDDFAL3-2-5-09  
 Client ID: ST020509MO      Analyte:      FAL ID: 1613 CS0 (080827J)

	Typ	Name	Amount	Resp	RA	RT	RF	RRF
1	Unk	2,3,7,8-TCDD	0.25	8.65e+04	0.75 y	27:27	-	1.06 y
2	Unk	1,2,3,7,8-PeCDD	1.25	2.81e+05	1.77 y	33:15	-	0.897 y
3	Unk	1,2,3,4,7,8-HxCDD	1.25	3.23e+05	1.26 y	38:36	-	1.45 y
4	Unk	1,2,3,6,7,8-HxCDD	1.25	2.46e+05	1.25 y	38:47	-	1.04 y
5	Unk	1,2,3,7,8,9-HxCDD	1.25	2.59e+05	1.20 y	39:12	-	1.13 y
6	Unk	1,2,3,4,6,7,8-HpCDD	1.25	2.09e+05	1.08 y	44:10	-	1.03 y
7	Unk	OCDD	2.50	2.79e+05	0.94 y	49:43	-	0.974 y
8	Unk	2,3,7,8-TCDF	0.25	1.85e+05	0.68 y	26:42	-	1.51 y
9	Unk	1,2,3,7,8-PeCDF	1.25	5.04e+05	1.42 y	31:32	-	1.01 y
10	Unk	2,3,4,7,8-PeCDF	1.25	4.39e+05	1.36 y	32:50	-	0.839 y
11	Unk	1,2,3,4,7,8-HxCDF	1.25	3.92e+05	1.26 y	37:12	-	1.29 y
12	Unk	1,2,3,6,7,8-HxCDF	1.25	3.65e+05	1.29 y	37:25	-	1.18 y
13	Unk	2,3,4,6,7,8-HxCDF	1.25	4.14e+05	1.23 y	38:20	-	1.46 y
14	Unk	1,2,3,7,8,9-HxCDF	1.25	2.64e+05	1.24 y	39:47	-	1.17 y
15	Unk	1,2,3,4,6,7,8-HpCDF	1.25	2.71e+05	0.97 y	42:18	-	1.04 y
16	Unk	1,2,3,4,7,8,9-HpCDF	1.25	2.10e+05	0.96 y	45:06	-	0.975 y
17	Unk	OCDF	2.50	3.63e+05	0.85 y	50:06	-	0.868 y
18	IS/RT	13C-2,3,7,8-TCDD	100.00	3.25e+07	0.79 y	27:25	-	0.908 y
19	IS	13C-1,2,3,7,8-PeCDD	100.00	2.50e+07	1.76 y	33:13	-	0.700 y
20	IS	13C-1,2,3,4,7,8-HxCDD	100.00	1.79e+07	1.25 y	38:35	-	0.993 y
21	IS	13C-1,2,3,6,7,8-HxCDD	100.00	1.89e+07	1.24 y	38:45	-	1.05 y
22	IS	13C-1,2,3,4,6,7,8-HpCDD	100.00	1.62e+07	1.11 y	44:10	-	0.903 y
23	IS	13C-OCDD	200.00	2.29e+07	0.91 y	49:42	-	0.637 y
24	IS	13C-2,3,7,8-TCDF	100.00	4.90e+07	0.82 y	26:41	-	0.848 y
25	IS	13C-1,2,3,7,8-PeCDF	100.00	4.00e+07	1.68 y	31:30	-	0.694 y
26	IS	13C-2,3,4,7,8-PeCDF	100.00	4.19e+07	1.69 y	32:49	-	0.725 y
27	IS	13C-1,2,3,4,7,8-HxCDF	100.00	2.43e+07	0.53 y	37:11	-	1.35 y
28	IS	13C-1,2,3,6,7,8-HxCDF	100.00	2.46e+07	0.53 y	37:23	-	1.37 y
29	IS	13C-2,3,4,6,7,8-HxCDF	100.00	2.26e+07	0.53 y	38:20	-	1.26 y
30	IS	13C-1,2,3,7,8,9-HxCDF	100.00	1.81e+07	0.53 y	39:46	-	1.01 y
31	IS	13C-1,2,3,4,6,7,8-HpCDF	100.00	2.09e+07	0.43 y	42:16	-	1.16 y
32	IS	13C-1,2,3,4,7,8,9-HpCDF	100.00	1.72e+07	0.42 y	45:05	-	0.958 y
33	IS	13C-OCDF	200.00	3.34e+07	0.91 y	50:05	-	0.930 y
34	C/Up	37Cl-2,3,7,8-TCDD	0.25	9.03e+04		27:26	-	1.01 y
35	RS	13C-1,2,3,4-TCDD	100.00	3.58e+07	0.80 y	26:50	3.58e+05	- n
36	RS	13C-1,2,3,4-TCDF	100.00	5.77e+07	0.82 y	25:35	5.77e+05	- n
37	RS/RT	13C-1,2,3,7,8,9-HxCDD	100.00	1.80e+07	1.25 y	39:11	1.80e+05	- n
38	Tot	Total Tetra-Dioxins	0.00	-	- n	-	-	1.06 y
39	Tot	Total Penta-Dioxins	0.00	-	- n	-	-	0.897 y
40	Tot	Total Hexa-Dioxins	0.00	-	- n	-	-	1.20 y
41	Tot	Total Hepta-Dioxins	0.00	-	- n	-	-	1.03 y
42	Tot	Total Tetra-Furans	0.00	-	- n	-	-	1.51 y
43	Tot	1st Fn. Tot Penta-Furans	0.00	-	- n	-	-	0.921 y
44	Tot	Total Penta-Furans	0.00	-	- n	-	-	0.921 y
45	Tot	Total Hexa-Furans	0.00	-	- n	-	-	1.28 y
46	Tot	Total Hepta-Furans	0.00	-	- n	-	-	1.01 y

Analyst: E

Date: 2/6/09

Run #2      Filename 05FEB09M      S: 3      Acquired: 5-FEB-09 17:10:30      Cal: PCDDFAL3-2-5-09  
 Client ID: ST020509M1      Analyte: FAL ID: 1613 CS1 (080827J)

	Typ	Name	Amount	Resp	RA	RT	RF	RRF
1	Unk	2,3,7,8-TCDD	0.50	1.95e+05	0.81 y	27:28	-	1.08 y
2	Unk	1,2,3,7,8-PeCDD	2.50	6.46e+05	1.64 y	33:15	-	0.903 y
3	Unk	1,2,3,4,7,8-HxCDD	2.50	7.35e+05	1.36 y	38:36	-	1.42 y
4	Unk	1,2,3,6,7,8-HxCDD	2.50	5.52e+05	1.25 y	38:46	-	1.04 y
5	Unk	1,2,3,7,8,9-HxCDD	2.50	5.65e+05	1.33 y	39:12	-	1.07 y
6	Unk	1,2,3,4,6,7,8-HpCDD	2.50	4.40e+05	1.10 y	44:11	-	1.01 y
7	Unk	OCDD	5.00	5.77e+05	0.91 y	49:43	-	0.939 y
8	Unk	2,3,7,8-TCDF	0.50	3.89e+05	0.68 y	26:42	-	1.43 y
9	Unk	1,2,3,7,8-PeCDF	2.50	1.15e+06	1.44 y	31:31	-	0.974 y
10	Unk	2,3,4,7,8-PeCDF	2.50	9.17e+05	1.38 y	32:50	-	0.757 y
11	Unk	1,2,3,4,7,8-HxCDF	2.50	8.71e+05	1.25 y	37:13	-	1.29 y
12	Unk	1,2,3,6,7,8-HxCDF	2.50	8.49e+05	1.28 y	37:25	-	1.22 y
13	Unk	2,3,4,6,7,8-HxCDF	2.50	9.66e+05	1.24 y	38:21	-	1.50 y
14	Unk	1,2,3,7,8,9-HxCDF	2.50	5.77e+05	1.28 y	39:47	-	1.15 y
15	Unk	1,2,3,4,6,7,8-HpCDF	2.50	6.17e+05	1.04 y	42:17	-	1.07 y
16	Unk	1,2,3,4,7,8,9-HpCDF	2.50	4.70e+05	1.00 y	45:06	-	0.999 y
17	Unk	OCDF	5.00	6.51e+05	0.91 y	50:07	-	0.710 y
18	IS/RT	13C-2,3,7,8-TCDD	100.00	3.62e+07	0.79 y	27:26	-	0.909 y
19	IS	13C-1,2,3,7,8-PeCDD	100.00	2.86e+07	1.76 y	33:14	-	0.718 y
20	IS	13C-1,2,3,4,7,8-HxCDD	100.00	2.08e+07	1.26 y	38:35	-	1.03 y
21	IS	13C-1,2,3,6,7,8-HxCDD	100.00	2.13e+07	1.26 y	38:45	-	1.06 y
22	IS	13C-1,2,3,4,6,7,8-HpCDD	100.00	1.75e+07	1.11 y	44:09	-	0.867 y
23	IS	13C-OCDD	200.00	2.46e+07	0.90 y	49:42	-	0.611 y
24	IS	13C-2,3,7,8-TCDF	100.00	5.45e+07	0.82 y	26:41	-	0.869 y
25	IS	13C-1,2,3,7,8-PeCDF	100.00	4.74e+07	1.68 y	31:31	-	0.755 y
26	IS	13C-2,3,4,7,8-PeCDF	100.00	4.84e+07	1.67 y	32:50	-	0.772 y
27	IS	13C-1,2,3,4,7,8-HxCDF	100.00	2.71e+07	0.53 y	37:12	-	1.35 y
28	IS	13C-1,2,3,6,7,8-HxCDF	100.00	2.79e+07	0.53 y	37:23	-	1.39 y
29	IS	13C-2,3,4,6,7,8-HxCDF	100.00	2.58e+07	0.54 y	38:19	-	1.28 y
30	IS	13C-1,2,3,7,8,9-HxCDF	100.00	2.00e+07	0.53 y	39:46	-	0.995 y
31	IS	13C-1,2,3,4,6,7,8-HpCDF	100.00	2.32e+07	0.43 y	42:15	-	1.15 y
32	IS	13C-1,2,3,4,7,8,9-HpCDF	100.00	1.88e+07	0.43 y	45:05	-	0.935 y
33	IS	13C-OCDF	200.00	3.67e+07	0.91 y	50:05	-	0.910 y
34	C/Up	37CL-2,3,7,8-TCDD	0.50	1.60e+05		27:28	-	0.805 y
35	RS	13C-1,2,3,4-TCDD	100.00	3.98e+07	0.79 y	26:51	3.98e+05	- n
36	RS	13C-1,2,3,4-TCDF	100.00	6.27e+07	0.82 y	25:36	6.27e+05	- n
37	RS/RT	13C-1,2,3,7,8,9-HxCDD	100.00	2.01e+07	1.25 y	39:11	2.01e+05	- n
38	Tot	Total Tetra-Dioxins	0.00	-	- n	-	-	1.08 y
39	Tot	Total Penta-Dioxins	0.00	-	- n	-	-	0.903 y
40	Tot	Total Hexa-Dioxins	0.00	-	- n	-	-	1.17 y
41	Tot	Total Hepta-Dioxins	0.00	-	- n	-	-	1.01 y
42	Tot	Total Tetra-Furans	0.00	-	- n	-	-	1.43 y
43	Tot	1st Fn. Tot Penta-Furans	0.00	-	- n	-	-	0.864 y
44	Tot	Total Penta-Furans	0.00	-	- n	-	-	0.864 y
45	Tot	Total Hexa-Furans	0.00	-	- n	-	-	1.29 y
46	Tot	Total Hepta-Furans	0.00	-	- n	-	-	1.04 y

Analyst: 

Date: 2/4/09

Run #3      Filename 05FEB09M      S: 4      Acquired: 5-FEB-09 18:05:53      Cal: PCDDFAL3-2-5-09  
 Client ID: ST020509M2      Analyte:      FAL ID: 1613 CS2 (080827J)

Typ	Name	Amount	Resp	RA	RT	RF	RRF
1	Unk                  2,3,7,8-TCDD	2.00	6.37e+05	0.71 y	27:26	-	0.980 y
2	Unk                  1,2,3,7,8-PeCDD	10.00	2.12e+06	1.66 y	33:16	-	0.870 y
3	Unk                  1,2,3,4,7,8-HxCDD	10.00	2.53e+06	1.30 y	38:36	-	1.40 y
4	Unk                  1,2,3,6,7,8-HxCDD	10.00	1.95e+06	1.32 y	38:47	-	1.00 y
5	Unk                  1,2,3,7,8,9-HxCDD	10.00	1.91e+06	1.29 y	39:12	-	1.02 y
6	Unk                  1,2,3,4,6,7,8-HpCDD	10.00	1.50e+06	1.04 y	44:11	-	0.962 y
7	Unk                  OCDD	20.00	1.95e+06	0.91 y	49:44	-	0.905 y
8	Unk                  2,3,7,8-TCDF	2.00	1.25e+06	0.71 y	26:42	-	1.29 y
9	Unk                  1,2,3,7,8-PeCDF	10.00	4.09e+06	1.47 y	31:32	-	1.01 y
10	Unk                  2,3,4,7,8-PeCDF	10.00	3.28e+06	1.45 y	32:50	-	0.797 y
11	Unk                  1,2,3,4,7,8-HxCDF	10.00	3.13e+06	1.26 y	37:13	-	1.30 y
12	Unk                  1,2,3,6,7,8-HxCDF	10.00	3.04e+06	1.22 y	37:25	-	1.21 y
13	Unk                  2,3,4,6,7,8-HxCDF	10.00	3.41e+06	1.26 y	38:21	-	1.46 y
14	Unk                  1,2,3,7,8,9-HxCDF	10.00	2.02e+06	1.27 y	39:47	-	1.13 y
15	Unk                  1,2,3,4,6,7,8-HpCDF	10.00	2.17e+06	1.03 y	42:17	-	1.05 y
16	Unk                  1,2,3,4,7,8,9-HpCDF	10.00	1.62e+06	1.02 y	45:06	-	0.966 y
17	Unk                  OCDF	20.00	2.46e+06	0.91 y	50:05	-	0.757 y
18	IS/RT                13C-2,3,7,8-TCDD	100.00	3.25e+07	0.78 y	27:25	-	0.940 y
19	IS                  13C-1,2,3,7,8-PeCDD	100.00	2.44e+07	1.78 y	33:14	-	0.706 y
20	IS                  13C-1,2,3,4,7,8-HxCDD	100.00	1.81e+07	1.28 y	38:35	-	1.03 y
21	IS                  13C-1,2,3,6,7,8-HxCDD	100.00	1.95e+07	1.25 y	38:45	-	1.11 y
22	IS                  13C-1,2,3,4,6,7,8-HpCDD	100.00	1.56e+07	1.12 y	44:10	-	0.888 y
23	IS                  13C-OCDD	200.00	2.16e+07	0.90 y	49:42	-	0.615 y
24	IS                  13C-2,3,7,8-TCDF	100.00	4.86e+07	0.82 y	26:41	-	0.855 y
25	IS                  13C-1,2,3,7,8-PeCDF	100.00	4.05e+07	1.70 y	31:30	-	0.713 y
26	IS                  13C-2,3,4,7,8-PeCDF	100.00	4.12e+07	1.69 y	32:49	-	0.724 y
27	IS                  13C-1,2,3,4,7,8-HxCDF	100.00	2.40e+07	0.53 y	37:11	-	1.37 y
28	IS                  13C-1,2,3,6,7,8-HxCDF	100.00	2.51e+07	0.53 y	37:24	-	1.43 y
29	IS                  13C-2,3,4,6,7,8-HxCDF	100.00	2.33e+07	0.53 y	38:19	-	1.33 y
30	IS                  13C-1,2,3,7,8,9-HxCDF	100.00	1.79e+07	0.54 y	39:45	-	1.02 y
31	IS                  13C-1,2,3,4,6,7,8-HpCDF	100.00	2.07e+07	0.42 y	42:16	-	1.18 y
32	IS                  13C-1,2,3,4,7,8,9-HpCDF	100.00	1.68e+07	0.42 y	45:05	-	0.957 y
33	IS                  13C-OCDF	200.00	3.25e+07	0.91 y	50:05	-	0.928 y
34	C/Up               37Cl-2,3,7,8-TCDD	2.00	5.74e+05		27:26	-	0.831 y
35	RS                  13C-1,2,3,4-TCDD	100.00	3.46e+07	0.79 y	26:51	3.46e+05	- n
36	RS                  13C-1,2,3,4-TCDF	100.00	5.68e+07	0.82 y	25:35	5.68e+05	- n
37	RS/RT               13C-1,2,3,7,8,9-HxCDD	100.00	1.75e+07	1.26 y	39:11	1.75e+05	- n
38	Tot                Total Tetra-Dioxins	0.00	-	- n	-	-	0.980 y
39	Tot                Total Penta-Dioxins	0.00	-	- n	-	-	0.870 y
40	Tot                Total Hexa-Dioxins	0.00	-	- n	-	-	1.13 y
41	Tot                Total Hepta-Dioxins	0.00	-	- n	-	-	0.962 y
42	Tot                Total Tetra-Furans	0.00	-	- n	-	-	1.29 y
43	Tot                1st Fn. Tot Penta-Furans	0.00	-	- n	-	-	0.902 y
44	Tot                Total Penta-Furans	0.00	-	- n	-	-	0.902 y
45	Tot                Total Hexa-Furans	0.00	-	- n	-	-	1.29 y
46	Tot                Total Hepta-Furans	0.00	-	- n	-	-	1.01 y

Analyst: SC

Date: 2/6/09

Run #4      Filename 05FEB09M      S: 1      Acquired: 5-FEB-09 15:20:00      Cal: PCDDFAL3-2-5-09  
 Client ID: ST020509M3      Analyte:      FAL ID: 1613 cs3 (080827J)

	Typ	Name	Amount	Resp	RA	RT	RF	RRF
1	Unk	2,3,7,8-TCDD	10.00	3.52e+06	0.76 y	27:26	-	0.982 y
2	Unk	1,2,3,7,8-PeCDD	50.00	1.28e+07	1.63 y	33:15	-	0.877 y
3	Unk	1,2,3,4,7,8-HxCDD	50.00	1.47e+07	1.32 y	38:37	-	1.38 y
4	Unk	1,2,3,6,7,8-HxCDD	50.00	1.06e+07	1.32 y	38:46	-	0.967 y
5	Unk	1,2,3,7,8,9-HxCDD	50.00	1.10e+07	1.30 y	39:12	-	1.02 y
6	Unk	1,2,3,4,6,7,8-HpCDD	50.00	8.87e+06	1.03 y	44:11	-	0.925 y
7	Unk	OCDD	100.00	1.25e+07	0.91 y	49:43	-	0.899 y
8	Unk	2,3,7,8-TCDF	10.00	6.70e+06	0.72 y	26:41	-	1.24 y
9	Unk	1,2,3,7,8-PeCDF	50.00	2.41e+07	1.53 y	31:31	-	1.03 y
10	Unk	2,3,4,7,8-PeCDF	50.00	1.86e+07	1.51 y	32:50	-	0.773 y
11	Unk	1,2,3,4,7,8-HxCDF	50.00	1.88e+07	1.26 y	37:13	-	1.33 y
12	Unk	1,2,3,6,7,8-HxCDF	50.00	1.79e+07	1.25 y	37:25	-	1.23 y
13	Unk	2,3,4,6,7,8-HxCDF	50.00	1.99e+07	1.26 y	38:21	-	1.51 y
14	Unk	1,2,3,7,8,9-HxCDF	50.00	1.18e+07	1.26 y	39:47	-	1.14 y
15	Unk	1,2,3,4,6,7,8-HpCDF	50.00	1.28e+07	1.03 y	42:17	-	1.07 y
16	Unk	1,2,3,4,7,8,9-HpCDF	50.00	1.01e+07	1.03 y	45:06	-	0.999 y
17	Unk	OCDF	100.00	1.40e+07	0.90 y	50:06	-	0.703 y
18	IS/RT	13C-2,3,7,8-TCDD	100.00	3.58e+07	0.78 y	27:25	-	0.914 y
19	IS	13C-1,2,3,7,8-PeCDD	100.00	2.93e+07	1.77 y	33:14	-	0.747 y
20	IS	13C-1,2,3,4,7,8-HxCDD	100.00	2.13e+07	1.27 y	38:35	-	1.03 y
21	IS	13C-1,2,3,6,7,8-HxCDD	100.00	2.20e+07	1.26 y	38:45	-	1.07 y
22	IS	13C-1,2,3,4,6,7,8-HpCDD	100.00	1.92e+07	1.10 y	44:10	-	0.928 y
23	IS	13C-OCDD	200.00	2.77e+07	0.92 y	49:42	-	0.671 y
24	IS	13C-2,3,7,8-TCDF	100.00	5.39e+07	0.81 y	26:40	-	0.864 y
25	IS	13C-1,2,3,7,8-PeCDF	100.00	4.67e+07	1.69 y	31:31	-	0.749 y
26	IS	13C-2,3,4,7,8-PeCDF	100.00	4.82e+07	1.70 y	32:50	-	0.772 y
27	IS	13C-1,2,3,4,7,8-HxCDF	100.00	2.81e+07	0.53 y	37:12	-	1.36 y
28	IS	13C-1,2,3,6,7,8-HxCDF	100.00	2.92e+07	0.53 y	37:23	-	1.41 y
29	IS	13C-2,3,4,6,7,8-HxCDF	100.00	2.64e+07	0.53 y	38:20	-	1.28 y
30	IS	13C-1,2,3,7,8,9-HxCDF	100.00	2.07e+07	0.54 y	39:46	-	1.00 y
31	IS	13C-1,2,3,4,6,7,8-HpCDF	100.00	2.39e+07	0.43 y	42:16	-	1.16 y
32	IS	13C-1,2,3,4,7,8,9-HpCDF	100.00	2.02e+07	0.42 y	45:05	-	0.977 y
33	IS	13C-OCDF	200.00	3.99e+07	0.91 y	50:05	-	0.967 y
34	C/Up	37Cl-2,3,7,8-TCDD	10.00	3.11e+06		27:26	-	0.794 y
35	RS	13C-1,2,3,4-TCDD	100.00	3.92e+07	0.79 y	26:51	3.92e+05	- n
36	RS	13C-1,2,3,4-TCDF	100.00	6.24e+07	0.82 y	25:35	6.24e+05	- n
37	RS/RT	13C-1,2,3,7,8,9-HxCDD	100.00	2.07e+07	1.27 y	39:12	2.07e+05	- n
38	Tot	Total Tetra-Dioxins	0.00	-	- n	-	-	0.982 y
39	Tot	Total Penta-Dioxins	0.00	-	- n	-	-	0.877 y
40	Tot	Total Hexa-Dioxins	0.00	-	- n	-	-	1.12 y
41	Tot	Total Hepta-Dioxins	0.00	-	- n	-	-	0.925 y
42	Tot	Total Tetra-Furans	0.00	-	- n	-	-	1.24 y
43	Tot	1st Fn. Tot Penta-Furans	0.00	-	- n	-	-	0.901 y
44	Tot	Total Penta-Furans	0.00	-	- n	-	-	0.901 y
45	Tot	Total Hexa-Furans	0.00	-	- n	-	-	1.31 y
46	Tot	Total Hepta-Furans	0.00	-	- n	-	-	1.04 y

Analyst:       

Date: 2/4/09

Run #5      Filename 05FEB09M      S: 5      Acquired: 5-FEB-09 19:01:11      Cal: PCDDFAL3-2-5-09  
 Client ID: ST020509M4      Analyte: FAL ID: 1613 CS4 (080827J)

	Typ	Name	Amount	Resp	RA	RT	RF	RRF
1	Unk	2,3,7,8-TCDD	40.00	1.56e+07	0.75 y	27:27	-	1.04 y
2	Unk	1,2,3,7,8-PeCDD	200.00	5.68e+07	1.63 y	33:16	-	0.917 y
3	Unk	1,2,3,4,7,8-HxCDD	200.00	6.57e+07	1.30 y	38:36	-	1.43 y
4	Unk	1,2,3,6,7,8-HxCDD	200.00	4.90e+07	1.31 y	38:46	-	1.03 y
5	Unk	1,2,3,7,8,9-HxCDD	200.00	4.72e+07	1.32 y	39:13	-	1.01 y
6	Unk	1,2,3,4,6,7,8-HpCDD	200.00	3.76e+07	1.04 y	44:11	-	0.979 y
7	Unk	OCDD	400.00	5.10e+07	0.90 y	49:44	-	0.954 y
8	Unk	2,3,7,8-TCDF	40.00	2.77e+07	0.73 y	26:42	-	1.25 y
9	Unk	1,2,3,7,8-PeCDF	200.00	1.07e+08	1.51 y	31:32	-	1.08 y
10	Unk	2,3,4,7,8-PeCDF	200.00	8.40e+07	1.53 y	32:51	-	0.826 y
11	Unk	1,2,3,4,7,8-HxCDF	200.00	8.43e+07	1.25 y	37:13	-	1.44 y
12	Unk	1,2,3,6,7,8-HxCDF	200.00	8.18e+07	1.25 y	37:25	-	1.31 y
13	Unk	2,3,4,6,7,8-HxCDF	200.00	9.09e+07	1.25 y	38:21	-	1.60 y
14	Unk	1,2,3,7,8,9-HxCDF	200.00	5.40e+07	1.25 y	39:47	-	1.24 y
15	Unk	1,2,3,4,6,7,8-HpCDF	200.00	5.85e+07	1.03 y	42:17	-	1.16 y
16	Unk	1,2,3,4,7,8,9-HpCDF	200.00	4.42e+07	1.05 y	45:06	-	1.08 y
17	Unk	OCDF	400.00	6.13e+07	0.91 y	50:06	-	0.755 y
18	IS/RT	13C-2,3,7,8-TCDD	100.00	3.76e+07	0.79 y	27:26	-	0.965 y
19	IS	13C-1,2,3,7,8-PeCDD	100.00	3.10e+07	1.77 y	33:14	-	0.795 y
20	IS	13C-1,2,3,4,7,8-HxCDD	100.00	2.29e+07	1.30 y	38:35	-	1.05 y
21	IS	13C-1,2,3,6,7,8-HxCDD	100.00	2.38e+07	1.24 y	38:45	-	1.09 y
22	IS	13C-1,2,3,4,6,7,8-HpCDD	100.00	1.92e+07	1.10 y	44:10	-	0.878 y
23	IS	13C-OCDD	200.00	2.67e+07	0.91 y	49:42	-	0.611 y
24	IS	13C-2,3,7,8-TCDF	100.00	5.54e+07	0.82 y	26:41	-	0.861 y
25	IS	13C-1,2,3,7,8-PeCDF	100.00	4.93e+07	1.71 y	31:30	-	0.765 y
26	IS	13C-2,3,4,7,8-PeCDF	100.00	5.09e+07	1.71 y	32:49	-	0.790 y
27	IS	13C-1,2,3,4,7,8-HxCDF	100.00	2.93e+07	0.53 y	37:12	-	1.34 y
28	IS	13C-1,2,3,6,7,8-HxCDF	100.00	3.12e+07	0.54 y	37:24	-	1.43 y
29	IS	13C-2,3,4,6,7,8-HxCDF	100.00	2.83e+07	0.53 y	38:19	-	1.30 y
30	IS	13C-1,2,3,7,8,9-HxCDF	100.00	2.17e+07	0.53 y	39:46	-	0.994 y
31	IS	13C-1,2,3,4,6,7,8-HpCDF	100.00	2.53e+07	0.43 y	42:16	-	1.16 y
32	IS	13C-1,2,3,4,7,8,9-HpCDF	100.00	2.04e+07	0.43 y	45:05	-	0.934 y
33	IS	13C-OCDF	200.00	4.06e+07	0.92 y	50:05	-	0.927 y
34	C/Up	37Cl-2,3,7,8-TCDD	40.00	1.37e+07		27:27	-	0.880 y
35	RS	13C-1,2,3,4-TCDD	100.00	3.90e+07	0.79 y	26:50	3.90e+05	- n
36	RS	13C-1,2,3,4-TCDF	100.00	6.44e+07	0.82 y	25:36	6.44e+05	- n
37	RS/RT	13C-1,2,3,7,8,9-HxCDD	100.00	2.19e+07	1.25 y	39:11	2.19e+05	- n
38	Tot	Total Tetra-Dioxins	0.00	-	- n	-	-	1.04 y
39	Tot	Total Penta-Dioxins	0.00	-	- n	-	-	0.917 y
40	Tot	Total Hexa-Dioxins	0.00	-	- n	-	-	1.15 y
41	Tot	Total Hepta-Dioxins	0.00	-	- n	-	-	0.979 y
42	Tot	Total Tetra-Furans	0.00	-	- n	-	-	1.25 y
43	Tot	1st Fn. Tot Penta-Furans	0.00	-	- n	-	-	0.953 y
44	Tot	Total Penta-Furans	0.00	-	- n	-	-	0.953 y
45	Tot	Total Hexa-Furans	0.00	-	- n	-	-	1.41 y
46	Tot	Total Hepta-Furans	0.00	-	- n	-	-	1.12 y

Analyst: E

Date: 01/16/09


Run #6      Filename 05FEB09M      S: 6      Acquired: 5-FEB-09 19:56:22      Cal: PCDDFAL3-2-5-09  
 Client ID: ST020509M5      Analyte:      FAL ID: 1613 CS5 (080827J)

	Typ	Name	Amount	Resp	RA	RT	RF	RRF
1	Unk	2,3,7,8-TCDD	200.00	7.46e+07	0.76 y	27:27	-	1.10 y
2	Unk	1,2,3,7,8-PeCDD	1000.00	2.79e+08	1.61 y	33:15	-	0.944 y
3	Unk	1,2,3,4,7,8-HxCDD	1000.00	3.45e+08	1.30 y	38:36	-	1.49 y
4	Unk	1,2,3,6,7,8-HxCDD	1000.00	2.55e+08	1.30 y	38:46	-	1.04 y
5	Unk	1,2,3,7,8,9-HxCDD	1000.00	2.52e+08	1.30 y	39:13	-	1.06 y
6	Unk	1,2,3,4,6,7,8-HpCDD	1000.00	2.04e+08	1.04 y	44:11	-	1.01 y
7	Unk	OCDD	2000.00	2.89e+08	0.90 y	49:45	-	0.969 y
8	Unk	2,3,7,8-TCDF	200.00	1.28e+08	0.73 y	26:41	-	1.31 y
9	Unk	1,2,3,7,8-PeCDF	1000.00	5.23e+08	1.52 y	31:32	-	1.08 y
10	Unk	2,3,4,7,8-PeCDF	1000.00	4.14e+08	1.51 y	32:51	-	0.858 y
11	Unk	1,2,3,4,7,8-HxCDF	1000.00	4.25e+08	1.26 y	37:13	-	1.47 y
12	Unk	1,2,3,6,7,8-HxCDF	1000.00	4.16e+08	1.24 y	37:25	-	1.35 y
13	Unk	2,3,4,6,7,8-HxCDF	1000.00	4.75e+08	1.24 y	38:21	-	1.65 y
14	Unk	1,2,3,7,8,9-HxCDF	1000.00	2.93e+08	1.26 y	39:47	-	1.28 y
15	Unk	1,2,3,4,6,7,8-HpCDF	1000.00	3.07e+08	1.03 y	42:17	-	1.21 y
16	Unk	1,2,3,4,7,8,9-HpCDF	1000.00	2.46e+08	1.04 y	45:06	-	1.12 y
17	Unk	OCDF	2000.00	3.60e+08	0.91 y	50:07	-	0.786 y
18	IS/RT	13C-2,3,7,8-TCDD	100.00	3.39e+07	0.79 y	27:25	-	0.984 y
19	IS	13C-1,2,3,7,8-PeCDD	100.00	2.96e+07	1.76 y	33:15	-	0.857 y
20	IS	13C-1,2,3,4,7,8-HxCDD	100.00	2.31e+07	1.26 y	38:35	-	1.03 y
21	IS	13C-1,2,3,6,7,8-HxCDD	100.00	2.46e+07	1.25 y	38:45	-	1.10 y
22	IS	13C-1,2,3,4,6,7,8-HpCDD	100.00	2.01e+07	1.10 y	44:10	-	0.898 y
23	IS	13C-OCDD	200.00	2.99e+07	0.91 y	49:43	-	0.667 y
24	IS	13C-2,3,7,8-TCDF	100.00	4.87e+07	0.82 y	26:39	-	0.880 y
25	IS	13C-1,2,3,7,8-PeCDF	100.00	4.84e+07	1.68 y	31:31	-	0.874 y
26	IS	13C-2,3,4,7,8-PeCDF	100.00	4.83e+07	1.68 y	32:50	-	0.871 y
27	IS	13C-1,2,3,4,7,8-HxCDF	100.00	2.88e+07	0.53 y	37:12	-	1.29 y
28	IS	13C-1,2,3,6,7,8-HxCDF	100.00	3.07e+07	0.54 y	37:23	-	1.37 y
29	IS	13C-2,3,4,6,7,8-HxCDF	100.00	2.88e+07	0.53 y	38:20	-	1.28 y
30	IS	13C-1,2,3,7,8,9-HxCDF	100.00	2.28e+07	0.54 y	39:46	-	1.02 y
31	IS	13C-1,2,3,4,6,7,8-HpCDF	100.00	2.54e+07	0.43 y	42:16	-	1.13 y
32	IS	13C-1,2,3,4,7,8,9-HpCDF	100.00	2.19e+07	0.43 y	45:05	-	0.980 y
33	IS	13C-OCDF	200.00	4.58e+07	0.92 y	50:06	-	1.02 y
34	C/Up	37Cl-2,3,7,8-TCDD	200.00	6.61e+07		27:27	-	0.958 y
35	RS	13C-1,2,3,4-TCDD	100.00	3.45e+07	0.80 y	26:50	3.45e+05	- n
36	RS	13C-1,2,3,4-TCDF	100.00	5.54e+07	0.83 y	25:35	5.54e+05	- n
37	RS/RT	13C-1,2,3,7,8,9-HxCDD	100.00	2.24e+07	1.24 y	39:12	2.24e+05	- n
38	Tot	Total Tetra-Dioxins	0.00	-	- n	-	-	1.10 y
39	Tot	Total Penta-Dioxins	0.00	-	- n	-	-	0.944 y
40	Tot	Total Hexa-Dioxins	0.00	-	- n	-	-	1.19 y
41	Tot	Total Hepta-Dioxins	0.00	-	- n	-	-	1.01 y
42	Tot	Total Tetra-Furans	0.00	-	- n	-	-	1.31 y
43	Tot	1st Fn. Tot Penta-Furans	0.00	-	- n	-	-	0.969 y
44	Tot	Total Penta-Furans	0.00	-	- n	-	-	0.969 y
45	Tot	Total Hexa-Furans	0.00	-	- n	-	-	1.45 y
46	Tot	Total Hepta-Furans	0.00	-	- n	-	-	1.17 y


Analyst: E

Date: 2/6/09


Peak Locate Examination: 5-FEB-2009:15:18 File:05FEB09M  
 Experiment:PCDD Function:1 Reference:PFK




Peak Locate Examination: 5-FEB-2009:15:18 File:05FEE09M  
 Experiment:PCDD Function:2 Reference:PFK




Peak Locate Examination: 5-FEB-2009:15:18 File:05FEB09M  
 Experiment:PCDD Function:3 Reference:PFK




Peak Locate Examination: 5-FEB-2009:15:19 File:05FEB09M  
 Experiment:PCDD Function:4 Reference:PTK




Peak Locate Examination: 5-FEB-2009:15:19 File:05FEB09M  
Experiment:PCDD Function:5 Reference:PFK




Peak Locate Examination: 6-FEB-2009:10:59 File:05FEB09M\_RES\_CHECK  
Experiment:PCDD Function:1 Reference:PFK




Peak Locate Examination: 6-FEB-2009:11:00 File:05FEB09M\_RES\_CHECK  
Experiment:PCDD Function:2 Reference:PFK




Peak Locate Examination: 6-FEB-2009:11:00 File:05FEB00M\_RES\_CHECK  
Experiment:PCDD Function:3 Reference:PFK



Peak Locate Examination: 6-FEB-2009:11:01 File:05FEB09M\_RES\_CHECK  
Experiment:PCDD Function:4 Reference:PFK



Peak Locate Examination: 6-FEB-2009:11:01 File:05FEB09M\_RES\_CHECK  
 Experiment:PCDD Function:5 Reference:PPK



File:05FEB09M #1-390 Acq: 5-FEB-2009 15:20:00 GC EI+ Voltage SIR Autospec-Ultima

319,8965 Exp:PCDD

Sample Text:ST020509M3 File Text:Frontier Analytical Laboratory

27:19

27:26

28:24

4.3E5

3.9E5

3.5E5

3.0E5

2.6E5

2.2E5


1.7E5

1.3E5

8.7E4

4.3E4

0.0E0



File:05FEB09M #1-390 Acq: 5-FEB-2009 15:20:00 GC EI+ Voltage SIR Autospec-Ultima

333,9339 Exp:PCDD

Sample Text:ST020509M3 File Text:Frontier Analytical Laboratory

26:51

27:25

5.5E6

5.0E6

4.4E6

3.9E6

3.3E6


2.8E6

2.2E6





1.7E6

1.1E6




5.5E5




File:05FEB09M #1-391 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima  
319.8965 S:2 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-424 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima  
355.8546 S:2 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-464 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima  
389.8156 S:2 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-464 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima  
391.8127 S:2 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory






File:05FEB09M #1-464 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima  
403.8530 S:2 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-464 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima  
380.9760 S:2 F:3 Exp:PCDD  
Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory

File:05FEB09M #1-541 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima  
423.7767 S:2 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-347 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima  
457.7377 S:2 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory

A1.35E5




File:05FEB09M #1-347 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima  
469.7780 S:2 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory

A1.09E7



File:05FEB09M #1-347 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima  
471.7750 S:2 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory

H1.67E6  
A1.20E7







File:05FEB09M #1-347 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima  
454.9728 S:2 F:5 Exp:PCDD  
Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory


100%  
48:23 48:46 49:07 49:25 49:40



File:05FEB09M #1-391 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima  
303.9016 S:2 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-391 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima  
 339.8597 S:2 BSUB(0.0000,15,-3.0) PKD(5.5,3.0,0.10%,100.0,0.00%,F,F) Exp:PCDD  
 Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-391 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima  
 341.8568 S:2 BSUB(10000,15,-3.0) PKD(5.5,3.0,0.10%,100.0,0.00%,F,F) Exp:PCDD  
 Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-391 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima  
 330.9792 S:2 Exp:PCDD  
 Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-424 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima  
339.8397 S:2 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,R,F) Exp:PCDD  
Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory  
100 %

A2.96E5

A2.53E5

6.8E4

3.4E4

0.0E0



File:05FEB09M #1-424 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima  
341.8368 S:2 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,R,F) Exp:PCDD  
Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory  
100 %

A2.08E5

A1.86E5

4.8E4

2.4E4

0.0E0



File:05FEB09M #1-424 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima  
351.9000 S:2 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,R,F) Exp:PCDD  
Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory  
100 %


A2.51E7

A2.63E7

5.8E6

2.9E6

0.0E0



File:05FEB09M #1-424 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima  
353.8970 S:2 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,R,F) Exp:PCDD  
Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory  
100 %

H3.41E6

H3.44E6

3.4E6

1.7E6

0.0E0



File:05FEB09M #1-424 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima  
409.7974 S:2 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,R,F) Exp:PCDD  
Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory  
100 %

A4.67E3

A4.30E3

2.9E3

1.5E3

0.0E0



File:05FEB09M #1-464 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima  
373.8207 S:2 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory

100 %

A2.19E5


A2.29E5

A1.46E5

4.5E4

2.2E4

0.0E0



File:05FEB09M #1-464 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima  
375.8178 S:2 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory

100 %

A1.73E5


A1.85E5

A1.18E5

3.6E4

1.8E4

0.0E0



File:05FEB09M #1-464 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima  
383.8639 S:2 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory

100 %

A1.73E5

A1.85E5

A1.18E5

3.6E4

1.8E4

0.0E0



File:05FEB09M #1-464 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima  
445.7555 S:2 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory

100 %

A1.73E5


A1.85E5

A1.18E5

3.6E4

1.8E4

0.0E0



File:05FEB09M #1-541 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima  
407.7818 S:2 F:4 BSUB(10000,15,-3,0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory


100 % A1.33E5

A1.03E5

2.3E4

1.2E4

0.0E0



File:05FEB09M #1-541 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima  
417.8253 S:2 F:4 BSUB(10000,15,-3,0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory

100 % A1.38E5

A1.07E5

2.4E4

1.2E4

0.0E0



File:05FEB09M #1-541 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima  
419.8220 S:2 F:4 BSUB(10000,15,-3,0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory


100 % A6.25E6

A5.07B6

1.1E6

5.6E5

0.0E0



File:05FEB09M #1-541 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima  
429.7165 S:2 F:4 BSUB(10000,15,-3,0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory


100 % H2.59E6

A1.46E7

2.6E6

1.3E6

0.0E0



File:05FEB09M #1-541 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima  
429.7165 S:2 F:4 BSUB(10000,15,-3,0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory


100 % A3.15E3

A1.76E3

1.7E3

1.3E6

0.0E0



File:05FEB09M #1-347 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima  
441.7428 S:2 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory


A1.66E5

A1.97E5

2.4E4

1.2E4

0.0E0



File:05FEB09M #1-347 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima

443.7398 S:2 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD


Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory

A1.39E7

2.9E4

1.5E4

0.0E0



File:05FEB09M #1-347 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima

453.7831 S:2 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD


Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory

A1.39E7

2.3E6

1.1E6

0.0E0



File:05FEB09M #1-347 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima

513.6775 S:2 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD


Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory

H2.51E6  
A1.75E7

2.5E6

1.3E6

0.0E0



File:05FEB09M #1-347 Acq: 5-FEB-2009 16:15:11 GC EI+ Voltage SIR Autospec-Ultima

513.6775 S:2 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD


Sample Text:ST020509M0 File Text:Frontier Analytical Laboratory

A891.16  
A249.82  
A173.44  
A140.22  
A129.22  
A81.63


5.1E2

2.5E2


0.0E0




File:05FEB09M #1-390 Acq: 5-FEB-2009 17:10:30 GC EI+ Voltage SIR Autospec-Ultima  
319.8965 S:3 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M1 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-390 Acq: 5-FEB-2009 17:10:30 GC EI+ Voltage SIR Autospec-Ultima  
321.8936 S:3 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M1 File Text:Frontier Analytical Laboratory







File:05FEB09M #1-390 Acq: 5-FEB-2009 17:10:30 GC EI+ Voltage SIR Autospec-Ultima  
331.9368 S:3 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M1 File Text:Frontier Analytical Laboratory







File:05FEB09M #1-390 Acq: 5-FEB-2009 17:10:30 GC EI+ Voltage SIR Autospec-Ultima  
333.9339 S:3 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M1 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-425 Acq: 5-FEB-2009 17:10:30 GC EI+ Voltage SIR Autospec-Ultima  
355.8546 S:3 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M1 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-463 Acq: 5-FEB-2009 17:10:30 GC EI+ Voltage SIR Autospec-Ultima  
389.8156 S:3 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M1 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-541 Acq: 5-FEB-2009 17:10:30 GC EI+ Voltage SIR Autospec-Ultima  
423.7767 S:3 F:4 BSUB(10000, 15,-3.0) PKD(5,5,3,0 10%,100.0,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M1 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-348 Acq: 5-FEB-2009 17:10:30 GC EI+ Voltage SIR Autospec-Ultima  
457.7377 S:3 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M1 File Text:Frontier Analytical Laboratory  
A2.74E5



File:05FEB09M #1-348 Acq: 5-FEB-2009 17:10:30 GC EI+ Voltage SIR Autospec-Ultima  
459.7348 S:3 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M1 File Text:Frontier Analytical Laboratory

A3.03E5

4.1E4

2.0E4

0.0E0

50.00

51.00

52.00

49:00

50:00

51:00

52:00

0

50

100 %

A1.17E7

File:05FEB09M #1-348 Acq: 5-FEB-2009 17:10:30 GC EI+ Voltage SIR Autospec-Ultima  
469.7780 S:3 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M1 File Text:Frontier Analytical Laboratory

A1.17E7

1.6E6

7.9E5

0.0E0

50.00

51.00

52.00

49:00

50:00

51:00

52:00

0

50

100 %

H1.78E6

A1.29E7

File:05FEB09M #1-348 Acq: 5-FEB-2009 17:10:30 GC EI+ Voltage SIR Autospec-Ultima  
471.7750 S:3 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M1 File Text:Frontier Analytical Laboratory

A1.29E7

1.8E6

8.9E5

0.0E0

50.00

51.00

52.00

49:00

50:00

51:00

52:00

0

50

100 %

5.7E7

2.9E7




0.0E0

Time

Time




Time

File:05FEB09M #1-390 Acq: 5-FEB-2009 17:10:30 GC EI+ Voltage SIR Autospec-Ultima  
303.9016 S:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M1 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-390 Acq: 5-FEB-2009 17:10:30 GC EI+ Voltage SIR Autospec-Ultima  
339.8597 S:3 BSUB(10000,15,-3,0) PKD(5,5,3,0,10%,10,0,0.00%,F,F) Exp:PCDD

Sample Text:ST020509M1 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-425 Acq: 5-FEB-2009 17:10:30 GC EI+ Voltage SIR Autospec-Ultima  
339.8597 S:3 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M1 File Text:Frontier Analytical Laboratory

100 %  
50  
0

1.6E5

7.9E4

0.0E0



File:05FEB09M #1-425 Acq: 5-FEB-2009 17:10:30 GC EI+ Voltage SIR Autospec-Ultima  
341.8568 S:3 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M1 File Text:Frontier Analytical Laboratory

100 %  
50  
0

A6.80E5

1.1E5

5.4E4

0.0E0

File:05FEB09M #1-425 Acq: 5-FEB-2009 17:10:30 GC EI+ Voltage SIR Autospec-Ultima  
351.9000 S:3 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M1 File Text:Frontier Analytical Laboratory

100 %  
50  
0

A4.73E5

A3.85E5

7.0E6

3.5E6

0.0E0

File:05FEB09M #1-425 Acq: 5-FEB-2009 17:10:30 GC EI+ Voltage SIR Autospec-Ultima  
353.8970 S:3 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M1 File Text:Frontier Analytical Laboratory

100 %  
50  
0

A2.97E7

A3.03E7

H4.09E6

A1.81E7

4.1E6

2.0E6

0.0E0

File:05FEB09M #1-425 Acq: 5-FEB-2009 17:10:30 GC EI+ Voltage SIR Autospec-Ultima  
409.7974 S:3 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M1 File Text:Frontier Analytical Laboratory

100 %  
50  
0

H4.00E6

A1.77E7

2.7E3

1.4E3

0.0E0

File:05FEB09M #1-425 Acq: 5-FEB-2009 17:10:30 GC EI+ Voltage SIR Autospec-Ultima  
409.7974 S:3 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M1 File Text:Frontier Analytical Laboratory

100 %  
50  
0

A4.45E3

A4.77E3

A2.84E3

A3.25E3

A3.67E3

A3.84E3

A3.73E3

A3.29E3

A5.39E3

A4.86E3

A3.86E3

A2.58E3

A4.92E3

A8.26E3

A4.18E3

A3.88E3

2.7E3

1.4E3

0.0E0

File:05FEB09M #1-425 Acq: 5-FEB-2009 17:10:30 GC EI+ Voltage SIR Autospec-Ultima  
409.7974 S:3 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M1 File Text:Frontier Analytical Laboratory

100 %  
50  
0

30:00

31:00

32:00

33:00


34:00

35:00

Time


File:05FEB09M #1-463 Acq: 5-FEB-2009 17:10:30 GC EI+ Voltage SIR Autospec-Ultima  
373.8207 S:3 F:3 BSUB(10000,15,-3,0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M1 File Text:Frontier Analytical Laboratory

A4.84E5  
A5.35E5  
A3.24E5  
8.2E4  
4.1E4  
0.0E0




File:05FEB09M #1-463 Acq: 5-FEB-2009 17:10:30 GC EI+ Voltage SIR Autospec-Ultima  
383.8639 S:3 F:3 BSUB(10000,15,-3,0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M1 File Text:Frontier Analytical Laboratory

A9.38E6  
A9.03E6  
A6.96E6  
1.8E6  
8.8E5  
0.0E0



H3.36E6  
H3.12E6  
A1.68E7  
H2.28E6  
A1.31E7  
3.4E6  
1.7E6  
0.0E0



File:05FEB09M #1-541 Acq: 5-FEB-2009 17:10:30 GC EI+ Voltage SIR Autospec-Ultima  
407.7818 S:3 F:4 BSUB(0000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M1 File Text:Frontier Analytical Laboratory

100 % A3.15E5

A2.35E5

5.6E4

2.8E4

0.0E0

42:00 43:00 44:00 45:00 46:00 47:00 48:00 Time

File:05FEB09M #1-541 Acq: 5-FEB-2009 17:10:30 GC EI+ Voltage SIR Autospec-Ultima  
409.7788 S:3 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD

Sample Text:ST020509M1 File Text:Frontier Analytical Laboratory

100 % A3.02E5

A2.36E5

5.1E4

2.6E4

0.0E0

42:00 43:00 44:00 45:00 46:00 47:00 48:00 Time

File:05FEB09M #1-541 Acq: 5-FEB-2009 17:10:30 GC EI+ Voltage SIR Autospec-Ultima  
417.8253 S:3 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD

Sample Text:ST020509M1 File Text:Frontier Analytical Laboratory

100 % A6.93E6

A5.67E6

1.2E6

6.0E5

0.0E0

42:00 43:00 44:00 45:00 46:00 47:00 48:00 Time

File:05FEB09M #1-541 Acq: 5-FEB-2009 17:10:30 GC EI+ Voltage SIR Autospec-Ultima  
419.8220 S:3 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD

Sample Text:ST020509M1 File Text:Frontier Analytical Laboratory

100 % H2.75E6

A1.63E7

H2.16E6

A1.32E7

2.7E6

1.4E6

0.0E0

42:00 43:00 44:00 45:00 46:00 47:00 48:00 Time

File:05FEB09M #1-541 Acq: 5-FEB-2009 17:10:30 GC EI+ Voltage SIR Autospec-Ultima  
479.7165 S:3 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD

Sample Text:ST020509M1 File Text:Frontier Analytical Laboratory

100 % A5.59E3

A4.56E3

A4.97E3

A3.31E3

A2.62E3

A2.74E3

A2.69E3

A3.82E3

A1.21E3

A1.03E3

A1.98E3

A1.39E3

A1.79E3


1.5E3

7.6E2

0.0E0


42:00 43:00 44:00 45:00 46:00 47:00 48:00 Time

File:05FEB09M #1-348 Acq: 5-FEB-2009 17:10:30 GC EI+ Voltage SIR Autospec-Ultima  
441.7428 S:3 F:5 BSUB(10000,15,-3,0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M1 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-348 Acq: 5-FEB-2009 17:10:30 GC EI+ Voltage SIR Autospec-Ultima  
443.7398 S:3 F:5 BSUB(10000,15,-3,0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M1 File Text:Frontier Analytical Laboratory

A3.40E5




File:05FEB09M #1-348 Acq: 5-FEB-2009 17:10:30 GC EI+ Voltage SIR Autospec-Ultima  
453.7831 S:3 F:5 BSUB(10000,15,-3,0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M1 File Text:Frontier Analytical Laboratory

A1.75E7






File:05FEB09M #1-348 Acq: 5-FEB-2009 17:10:30 GC EI+ Voltage SIR Autospec-Ultima  
455.7801 S:3 F:5 BSUB(10000,15,-3,0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M1 File Text:Frontier Analytical Laboratory

H2.57E6  
A1.92E7




File:05FEB09M #1-390 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
319.8965 S:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-425 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
355.8546 S:4 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-425 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
357.8517 S:4 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-425 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
367.8949 S:4 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-425 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
369.8919 S:4 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-425 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
369.8919 S:4 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-425 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
366.9792 S:4 F:2 Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-464 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
389.8156 S:4 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-464 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
401.8559 S:4 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-464 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
403.8530 S:4 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-541 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
423.7767 S:4 F:4 BSUB(10000,15,-3) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-541 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
425.7737 S:4 F:4 BSUB(10000,15,-3) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-541 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
435.8169 S:4 F:4 BSUB(10000,15,-3) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-541 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
437.8140 S:4 F:4 BSUB(10000,15,-3) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-541 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
430.9728 S:4 F:4 Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-347 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
457.7377 S:4 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory

100 %  
A9.31E5




File:05FEB09M #1-347 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
469.7780 S:4 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory

A1.02E7




File:05FEB09M #1-347 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
454.9728 S:4 F:5 Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory


100 %  
48.19  
48.38




File:05FEB09M #1-390 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
303.9016 S:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-390 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
315.9419 S:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-390 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
317.9389 S:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-390 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima

339.8597 S:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD


Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-390 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima

341.8568 S:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD


Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-390 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima

409.7974 S:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD


Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-390 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima

330.9792 S:4 Exp:PCDD

Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-425 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
339.8597 S:4 F:2 BSUB(100000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory

A2.43E6

A1.94E6

5.5E5

2.7E5



File:05FEB09M #1-425 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
341.8568 S:4 F:2 BSUB(100000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory

A1.65E6

A1.34E6

3.6E5

1.8E5




File:05FEB09M #1-425 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
351.9000 S:4 F:2 BSUB(100000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory

A2.55E7

A2.59E7

5.8E6

2.9E6



File:05FEB09M #1-425 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
353.8970 S:4 F:2 BSUB(100000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory

H3.42E6

H3.42E6

3.4E6

1.7E6




File:05FEB09M #1-425 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
409.7974 S:4 F:2 BSUB(100000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory

A7.00E3

A7.00E3

3.6E3

1.7E3



File:05FEB09M #1-464 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
373.8207 S:4 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD

Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory

100 %

3.6E5

1.8E5

0.0E0

50 %

2.9E5

1.5E5

0.0E0

0 %

A1.13E6

3.6E5

1.8E5

0.0E0

File:05FEB09M #1-464 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
375.8178 S:4 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD

Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory

100 %

A1.90E6

3.6E5

1.8E5

50 %

A1.38E6

2.9E5

1.5E5

0.0E0

0 %

A1.51E6

2.9E5

1.5E5

0.0E0

File:05FEB09M #1-464 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
383.8639 S:4 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD

Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory

100 %

A8.74E6

2.9E5

1.5E5

0.0E0

0 %

A8.06E6

2.9E5

1.5E5

0.0E0

File:05FEB09M #1-464 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
385.8610 S:4 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD

Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory

100 %

H3.14E6

2.9E5

1.5E5

0.0E0

0 %

A1.64E7

2.9E5

1.5E5

0.0E0

File:05FEB09M #1-464 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
445.7555 S:4 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD

Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory

100 %

A2.17E3

2.0E3

1.0E3

0.0E0

50 %

A1.19E3

9.9E2

1.0E3

0.0E0

0 %

A1.81E3

9.9E2

1.0E3

0.0E0

0 %

A1.10E3

9.9E2

1.0E3

0.0E0

0 %

A1.39E3

9.9E2

1.0E3

0.0E0

0 %

A1.07E3

9.9E2

1.0E3

0.0E0

0 %

A1.17E3

9.9E2

1.0E3

0.0E0

0 %






A1.64E3

9.9E2


1.0E3

0.0E0


File:05FEB09M #1-541 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
409.7788 S:4 F:4 BSUB(10000, 15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-347 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
441.7428 S:4 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-347 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
443.7398 S:4 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-347 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
453.7831 S:4 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory

A1.55E7

H2.37E6

A1.70E7

2.2E6

1.1E6

0.0E0

1.9E5

9.3E4

8.1E4

1.6E5

0.0E0

Time

File:05FEB09M #1-347 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
455.7801 S:4 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory

A1.55E7

H2.37E6

A1.70E7

2.4E6

1.2E6

0.0E0

1.9E5

9.3E4

8.1E4

1.6E5

0.0E0

Time

File:05FEB09M #1-347 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
513.6775 S:4 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory

A1.53E3

H2.37E6

A1.70E7

1.0E3

5.2E2

0.0E0

1.9E5

9.3E4

8.1E4

1.6E5

0.0E0

Time

File:05FEB09M #1-347 Acq: 5-FEB-2009 18:05:53 GC EI+ Voltage SIR Autospec-Ultima  
513.6775 S:4 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M2 File Text:Frontier Analytical Laboratory

A1.53E3

H2.37E6

A1.70E7

1.0E3

5.2E2

0.0E0

1.9E5

9.3E4



8.1E4

1.6E5

0.0E0


Time

File:05FEB09M #1-390 Acq: 5-FEB-2009 15:20:00 GC EI+ Voltage SIR Autospec-Ultima  
319.8965 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M3 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-425 Acq: 5-FEB-2009 15:20:00 GC EI+ Voltage SIR Autospec-Ultima  
355.8546 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M3 File Text:Frontier Analytical Laboratory

100 %  
A.1.02E7




File:05FEB09M #1-425 Acq: 5-FEB-2009 15:20:00 GC EI+ Voltage SIR Autospec-Ultima  
357.8517 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M3 File Text:Frontier Analytical Laboratory

100 %  
A6.29E6




File:05FEB09M #1-425 Acq: 5-FEB-2009 15:20:00 GC EI+ Voltage SIR Autospec-Ultima  
367.8949 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M3 File Text:Frontier Analytical Laboratory

100 %  
A1.87E7



File:05FEB09M #1-425 Acq: 5-FEB-2009 15:20:00 GC EI+ Voltage SIR Autospec-Ultima  
369.8919 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M3 File Text:Frontier Analytical Laboratory

100 %  
H2.35E6  
A1.06E7




File:05FEB09M #1-425 Acq: 5-FEB-2009 15:20:00 GC EI+ Voltage SIR Autospec-Ultima  
366.9792 F:2 Exp:PCDD  
Sample Text:ST020509M3 File Text:Frontier Analytical Laboratory


100 %  
29:40 29:58 30:23 30:57 31:24 31:48 32:18 32:45 33:10 33:31 33:58 34:29 34:56 3.7E7

0 50 100 %  
30:00 31:00 32:00 33:00 34:00 35:00 Time  
1.8E7  
0.0E0


File:05FEB09M #1-463 Acq: 5-FEB-2009 15:20:00 GC EI+ Voltage SIR Autospec-Ultima  
389.8156 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M3 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-463 Acq: 5-FEB-2009 15:20:00 GC EI+ Voltage SIR Autospec-Ultima  
391.8127 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M3 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-463 Acq: 5-FEB-2009 15:20:00 GC EI+ Voltage SIR Autospec-Ultima  
401.8559 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M3 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-463 Acq: 5-FEB-2009 15:20:00 GC EI+ Voltage SIR Autospec-Ultima  
380.9760 F:3 Exp:PCDD  
Sample Text:ST020509M3 File Text:Frontier Analytical Laboratory

File:05FEB09M #1-541 Acq: 5-FEB-2009 15:20:00 GC EI+ Voltage SIR Autospec-Ultima  
423.7767 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M3 File Text:Frontier Analytical Laboratory

100 %  
A6.61E6

A6.32E6

A4.50E6

A4.37E6

1.3E6

6.3E5

0.0E0



File:05FEB09M #1-541 Acq: 5-FEB-2009 15:20:00 GC EI+ Voltage SIR Autospec-Ultima  
425.7737 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M3 File Text:Frontier Analytical Laboratory


100 %  
A6.32E6

A4.37E6

1.2E6

6.0E5

0.0E0



File:05FEB09M #1-541 Acq: 5-FEB-2009 15:20:00 GC EI+ Voltage SIR Autospec-Ultima  
435.8169 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M3 File Text:Frontier Analytical Laboratory

100 %  
A1.00E7

A1.00E7

1.9E6

9.6E5

0.0E0



File:05FEB09M #1-541 Acq: 5-FEB-2009 15:20:00 GC EI+ Voltage SIR Autospec-Ultima  
437.8140 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M3 File Text:Frontier Analytical Laboratory

100 %  
H1.78E6

H1.78E6

1.8E6

8.9E5

0.0E0



File:05FEB09M #1-541 Acq: 5-FEB-2009 15:20:00 GC EI+ Voltage SIR Autospec-Ultima  
430.9728 F:4 Exp:PCDD  
Sample Text:ST020509M3 File Text:Frontier Analytical Laboratory


100 %  
A9.14E6

A9.14E6





3.9E7

2.0E7




0.0E0




File:05FEB09M #1-348 Acq: 5-FEB-2009 15:20:00 GC EI+ Voltage SIR Autospec-Ultima  
457.7377 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M3 File Text:Frontier Analytical Laboratory  
100 %



File:05FEB09M #1-390 Acq: 5-FEB-2009 15:20:00 GC EI+ Voltage SIR Autospec-Ultima  
303.9016 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M3 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-390 Acq: 5-FEB-2009 15:20:00 GC EI+ Voltage SIR Autospec-Ultima  
339.8597 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,FF) Exp:PCDD  
Sample Text:ST020509M3 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-425 Acq: 5-FEB-2009 15:20:00 GC EI+ Voltage SIR Autospec-Ultima  
 339.8597 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
 Sample Text:ST020509M3 File Text:Frontier Analytical Laboratory

A1.46E7

A1.12E7

A1.13E7

3.4E6

1.7E6

0.0E0

File:05FEB09M #1-425 Acq: 5-FEB-2009 15:20:00 GC EI+ Voltage SIR Autospec-Ultima  
 341.8568 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
 Sample Text:ST020509M3 File Text:Frontier Analytical Laboratory

A9.54E6

A7.41E6

A7.45E6

2.2E6

1.1E6

0.0E0

File:05FEB09M #1-425 Acq: 5-FEB-2009 15:20:00 GC EI+ Voltage SIR Autospec-Ultima  
 351.9000 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
 Sample Text:ST020509M3 File Text:Frontier Analytical Laboratory

A2.94E7

A3.03E7

6.9E6

3.4E6

0.0E0

File:05FEB09M #1-425 Acq: 5-FEB-2009 15:20:00 GC EI+ Voltage SIR Autospec-Ultima  
 353.8970 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
 Sample Text:ST020509M3 File Text:Frontier Analytical Laboratory

H4.02E6

H4.04E6

4.0E6

2.0E6

0.0E0

File:05FEB09M #1-425 Acq: 5-FEB-2009 15:20:00 GC EI+ Voltage SIR Autospec-Ultima  
 409.7974 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
 Sample Text:ST020509M3 File Text:Frontier Analytical Laboratory

A7.60E3

A7.31E3

A5.41E3

3.3E3

2.0E6

0.0E0

File:05FEB09M #1-425 Acq: 5-FEB-2009 15:20:00 GC EI+ Voltage SIR Autospec-Ultima  
 409.7974 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
 Sample Text:ST020509M3 File Text:Frontier Analytical Laboratory

A9.13E3

A3.00E3

A3.67E3

1.6E3

0.0E0

0.0E0

A3.81E3

A5.03E3






A3.58E3

1.6E3


0.0E0

0.0E0

File:05FEB09M #1-463 Acq: 5-FEB-2009 15:20:00 GC EI+ Voltage SIR Autospec-Ultima  
 373.8207 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
 Sample Text:ST020509M3 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-541 Acq: 5-FEB-2009 15:20:00 GC EI+ Voltage SIR Autospec-Ultima  
407.7818 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M3 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-541 Acq: 5-FEB-2009 15:20:00 GC EI+ Voltage SIR Autospec-Ultima  
417.8253 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M3 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-541 Acq: 5-FEB-2009 15:20:00 GC EI+ Voltage SIR Autospec-Ultima  
419.8220 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M3 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-541 Acq: 5-FEB-2009 15:20:00 GC EI+ Voltage SIR Autospec-Ultima  
479.7165 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M3 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-348 Acq: 5-FEB-2009 15:20:00 GC EI+ Voltage SIR Autospec-Ultima  
441.7428 F,5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M3 File Text:Frontier Analytical Laboratory

1.0E6  
5.2E5  
0.0E0




A7.37E6



A1.91E7




H3.14E6  
A2.09E7



3.1E6  
1.6E6  
0.0E0







A815.90  
A352.65  
A440.41  
A795.64  
A873.52  
6.8E2  
3.4E2  
0.0E0




A776.94  
A228.50  
A273.21  
A655.44  
A343.30  
A443.00  
A307.75  
A240.61  
A611.84  
A725.63  
A440.41  
A815.90  
A352.65  
A795.64  
A873.52  
6.8E2  
3.4E2  
0.0E0


File:05FEB09M #1-390 Acq: 5-FEB-2009 19:01:11 GC EI+ Voltage SIR Autospec-Ultima  
319.8965 S:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M4 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-424 Acq: 5-FEB-2009 19:01:11 GC EI+ Voltage SIR Autospec-Ultima  
355.8546 S:5 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M4 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-424 Acq: 5-FEB-2009 19:01:11 GC EI+ Voltage SIR Autospec-Ultima  
357.8517 S:5 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M4 File Text:Frontier Analytical Laboratory







File:05FEB09M #1-424 Acq: 5-FEB-2009 19:01:11 GC EI+ Voltage SIR Autospec-Ultima  
367.8949 S:5 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M4 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-424 Acq: 5-FEB-2009 19:01:11 GC EI+ Voltage SIR Autospec-Ultima  
369.8919 S:5 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M4 File Text:Frontier Analytical Laboratory







File:05FEB09M #1-464 Acq: 5-FEB-2009 19:01:11 GC EI+ Voltage SIR Autospec-Ultima  
389.8156 S:5 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M4 File Text:Frontier Analytical Laboratory







File:05FEB09M #1-541 Acq: 5-FEB-2009 19:01:11 GC EI+ Voltage SIR Autospec-Ultima  
423.7767 S:5 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M4 File Text:Frontier Analytical Laboratory





A1.92E7



File:05FEB09M #1-347 Acq: 5-FEB-2009 19:01:11 GC EI+ Voltage SIR Autospec-Ultima  
457.7377 S:5 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M4 File Text:Frontier Analytical Laboratory  
100 %  
A2.42E7




File:05FEB09M #1-390 Acq: 5-FEB-2009 19:01:11 GC EI+ Voltage SIR Autospec-Ultima  
303.9016 S:5 BSUB(0000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M4 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-390 Acq: 5-FEB-2009 19:01:11 GC EI+ Voltage SIR Autospec-Ultima  
339.8597 S:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD

Sample Text:ST020509M4 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-390 Acq: 5-FEB-2009 19:01:11 GC EI+ Voltage SIR Autospec-Ultima  
341.8568 S:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD

Sample Text:ST020509M4 File Text:Frontier Analytical Laboratory







File:05FEB09M #1-390 Acq: 5-FEB-2009 19:01:11 GC EI+ Voltage SIR Autospec-Ultima  
409.7974 S:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD

Sample Text:ST020509M4 File Text:Frontier Analytical Laboratory








File:05FEB09M #1-424 Acq: 5-FEB-2009 19:01:11 GC EI+ Voltage SIR Autospec-Ultima  
339.8597 S:5 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M4 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-464 Acq: 5-FEB-2009 19:01:11 GC EI+ Voltage SIR Autospec-Ultima  
373.8207 S:5 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M4 File Text:Frontier Analytical Laboratory


100 %  
A4.67E7  
A5.04E7  
A3.00E7  
-9.9E6  
-5.0E6  
0.0E0




File:05FEB09M #1-541 Acq: 5-FEB-2009 19:01:11 GC EI+ Voltage SIR Autospec-Ultima  
407.7818 S:5 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M4 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-541 Acq: 5-FEB-2009 19:01:11 GC EI+ Voltage SIR Autospec-Ultima  
409.7788 S:5 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M4 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-541 Acq: 5-FEB-2009 19:01:11 GC EI+ Voltage SIR Autospec-Ultima  
417.8253 S:5 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M4 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-541 Acq: 5-FEB-2009 19:01:11 GC EI+ Voltage SIR Autospec-Ultima  
419.8220 S:5 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M4 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-541 Acq: 5-FEB-2009 19:01:11 GC EI+ Voltage SIR Autospec-Ultima  
479.7165 S:5 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M4 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-347 Acq: 5-FEB-2009 19:01:11 GC EI+ Voltage SIR Autospec-Ultima  
441.7428 S:5 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M4 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-347 Acq: 5-FEB-2009 19:01:11 GC EI+ Voltage SIR Autospec-Ultima  
443.7398 S:5 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M4 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-347 Acq: 5-FEB-2009 19:01:11 GC EI+ Voltage SIR Autospec-Ultima  
453.7831 S:5 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M4 File Text:Frontier Analytical Laboratory

A1.94E7

2.7E6

1.3E6

0.0E0




File:05FEB09M #1-347 Acq: 5-FEB-2009 19:01:11 GC EI+ Voltage SIR Autospec-Ultima  
513.6775 S:5 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M4 File Text:Frontier Analytical Laboratory

A1.41E3





1.1E3

5.7E2


0.0E0




File:05FEB09M #1-390 Acq: 5-FEB-2009 19:56:22 GC EI+ Voltage SIR Autospec-Ultima  
319.8965 S:6 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M5 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-425 Acq: 5-FEB-2009 19:56:22 GC EI+ Voltage SIR Autospec-Ultima  
355.8546 S:6 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509MS File Text:Frontier Analytical Laboratory




File:05FEB09M #1-425 Acq: 5-FEB-2009 19:56:22 GC EI+ Voltage SIR Autospec-Ultima  
367.8949 S:6 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509MS File Text:Frontier Analytical Laboratory




File:05FEB09M #1-425 Acq: 5-FEB-2009 19:56:22 GC EI+ Voltage SIR Autospec-Ultima  
369.8919 S:6 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509MS File Text:Frontier Analytical Laboratory




File:05FEB09M #1-425 Acq: 5-FEB-2009 19:56:22 GC EI+ Voltage SIR Autospec-Ultima  
366.9792 S:6 F:2 Exp:PCDD  
Sample Text:ST020509MS File Text:Frontier Analytical Laboratory




File:05FEB09M #1-425 Acq: 5-FEB-2009 19:56:22 GC EI+ Voltage SIR Autospec-Ultima  
366.9792 S:6 F:2 Exp:PCDD  
Sample Text:ST020509MS File Text:Frontier Analytical Laboratory




File:05FEB09M #1-463 Acq: 5-FEB-2009 19:56:22 GC EI+ Voltage SIR Autospec-Ultima  
389.8156 S:6 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509MS File Text:Frontier Analytical Laboratory






File:05FEB09M #1-463 Acq: 5-FEB-2009 19:56:22 GC EI+ Voltage SIR Autospec-Ultima  
401.8559 S:6 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509MS File Text:Frontier Analytical Laboratory



File:05FEB09M #1-463 Acq: 5-FEB-2009 19:56:22 GC EI+ Voltage SIR Autospec-Ultima  
380.9760 S:6 F:3 Exp:PCDD  
Sample Text:ST020509MS File Text:Frontier Analytical Laboratory







File:05FEB09M #1-541 Acq: 5-FEB-2009 19:56:22 GC EI+ Voltage SIR Autospec-Ultima  
423.7767 S:6 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M5 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-348 Acq: 5-FEB-2009 19:56:22 GC EI+ Voltage SIR Autospec-Ultima  
457.7377 S:6 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509MS File Text:Frontier Analytical Laboratory


100 %  
A1.37E8



File:05FEB09M #1-390 Acq: 5-FEB-2009 19:56:22 GC EI+ Voltage SIR Autospec-Ultima  
303.9016 S:6 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509MS File Text:Frontier Analytical Laboratory




File:05FEB09M #1-390 Acq: 5-FEB-2009 19:56:22 GC EI+ Voltage SIR Autospec-Ultima  
305.8987 S:6 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509MS File Text:Frontier Analytical Laboratory



File:05FEB09M #1-390 Acq: 5-FEB-2009 19:56:22 GC EI+ Voltage SIR Autospec-Ultima  
315.9419 S:6 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509MS File Text:Frontier Analytical Laboratory



File:05FEB09M #1-390 Acq: 5-FEB-2009 19:56:22 GC EI+ Voltage SIR Autospec-Ultima  
317.9389 S:6 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509MS File Text:Frontier Analytical Laboratory



File:05FEB09M #1-390 Acq: 5-FEB-2009 19:56:22 GC EI+ Voltage SIR Autospec-Ultima  
 339.8597 S:6 BSUB(0000,15,-3,0) PKD(5,5,3,0,10%,100,0,0.00%,F,F) Exp:PCDD  
 Sample Text:ST020509MS File Text:Frontier Analytical Laboratory

100 %  
 A4.87E3

1.4E3  
 1.1E3

8.3E2

5.5E2

2.8E2

0.0E0



File:05FEB09M #1-390 Acq: 5-FEB-2009 19:56:22 GC EI+ Voltage SIR Autospec-Ultima  
 341.8568 S:6 BSUB(10000,15,-3,0) PKD(5,5,3,0,10%,100,0,0.00%,F,F) Exp:PCDD  
 Sample Text:ST020509MS File Text:Frontier Analytical Laboratory

100 %  
 A8.81E3 A7.58E3 A5.37E3 A7.74E3 A7.47E3 A5.86E3 A7.98E3 A1.10E4 A1.04E4 A8.62E3 A1.27E4 A4.47E3 A8.64E3 A1.19E4 2.4E3 2.0E3 1.5E3 9.8E2 4.9E2 0.0E0


2.6E3

2.0E3


1.3E3

6.6E2

0.0E0




File:05FEB09M #1-390 Acq: 5-FEB-2009 19:56:22 GC EI+ Voltage SIR Autospec-Ultima  
 330.9792 S:6 Exp:PCDD  
 Sample Text:ST020509MS File Text:Frontier Analytical Laboratory




File:05FEB09M #1-425 Acq: 5-FEB-2009 19:56:22 GC El+ Voltage SIR Autospec-Ultima  
 339.8597 S:6 F:2 BSUB(10000,15,3.0) PKD(5,5,3,0,10%,100,0,0.0%,F,F) Exp:PCDD  
 Sample Text:ST020509M5 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-425 Acq: 5-FEB-2009 19:56:22 GC El+ Voltage SIR Autospec-Ultima  
 341.8568 S:6 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100,0,0.0%,F,F) Exp:PCDD  
 Sample Text:ST020509M5 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-425 Acq: 5-FEB-2009 19:56:22 GC El+ Voltage SIR Autospec-Ultima  
 351.9000 S:6 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100,0,0.0%,F,F) Exp:PCDD  
 Sample Text:ST020509M5 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-425 Acq: 5-FEB-2009 19:56:22 GC El+ Voltage SIR Autospec-Ultima  
 353.8970 S:6 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100,0,0.0%,F,F) Exp:PCDD  
 Sample Text:ST020509M5 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-463 Acq: 5-FEB-2009 19:56:22 GC El+ Voltage SIR Autospec-Ultima  
373.8207 S:6 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100,0,0.0%,F,F) Exp:PCDD

Sample Text:ST020509M5 File Text:Frontier Analytical Laboratory

A2.37E8

A2.63E8

A2.12E8

A1.63E8

5.2E7

2.6E7

0.0E0



File:05FEB09M #1-463 Acq: 5-FEB-2009 19:56:22 GC El+ Voltage SIR Autospec-Ultima  
383.8639 S:6 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100,0,0.0%,F,F) Exp:PCDD

Sample Text:ST020509M5 File Text:Frontier Analytical Laboratory

A1.88E8

A2.12E8

A1.30E8

4.1E7

2.1E7

0.0E0



File:05FEB09M #1-463 Acq: 5-FEB-2009 19:56:22 GC El+ Voltage SIR Autospec-Ultima  
385.8610 S:6 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100,0,0.0%,F,F) Exp:PCDD

Sample Text:ST020509M5 File Text:Frontier Analytical Laboratory

A1.08E7

A9.98E6

A8.01E6

2.1E6

1.0E6

0.0E0



File:05FEB09M #1-463 Acq: 5-FEB-2009 19:56:22 GC El+ Voltage SIR Autospec-Ultima  
445.7555 S:6 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100,0,0.0%,F,F) Exp:PCDD

Sample Text:ST020509M5 File Text:Frontier Analytical Laboratory

A2.89E3


A2.75E3

H3.87E6





1.6E3

8.2E2

0.0E0



File:05FEB09M #1-541 Acq: 5-FEB-2009 19:56:22 GC EI+ Voltage SIR Autospec-Ultima  
 407.7818 S:6 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
 Sample Text:ST020509M5 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-348 Acq: 5-FEB-2009 19:56:22 GC EI+ Voltage SIR Autospec-Ultima  
441.7428 S:6 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M5 File Text:Frontier Analytical Laboratory




File:05FEB09M #1-348 Acq: 5-FEB-2009 19:56:22 GC EI+ Voltage SIR Autospec-Ultima  
443.7398 S:6 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M5 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-348 Acq: 5-FEB-2009 19:56:22 GC EI+ Voltage SIR Autospec-Ultima  
453.7831 S:6 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M5 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-348 Acq: 5-FEB-2009 19:56:22 GC EI+ Voltage SIR Autospec-Ultima  
455.7801 S:6 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M5 File Text:Frontier Analytical Laboratory



File:05FEB09M #1-348 Acq: 5-FEB-2009 19:56:22 GC EI+ Voltage SIR Autospec-Ultima  
513.6775 S:6 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST020509M5 File Text:Frontier Analytical Laboratory



6DFA - FORM VI-HR CDD-1  
TCDF INITIAL CALIBRATION RESPONSE FACTOR SUMMARY  
HIGH RESOLUTION

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: SDG NO.: SDG E3YC9

GC COLUMN: DB225 ID: 0.25 (mm) INSTRUMENT ID: FAL1

INIT. CALIB. DATE(S): 13-JAN-09 13-JAN-09 13-JAN-09 13-JAN-09 13-JAN-09

INIT. CALIB. TIMES: 11:41:31 12:16:35 11:06:28 12:51:39 13:26:43

TARGET ANALYTES	RR/RRF					MEAN RR/RRF	%RSD	QC LIMITS
	CS1	CS2	CS3	CS4	CS5			
2,3,7,8-TCDF	1.12	0.97	1.06	1.00	0.98	1.03	6.35	+/-20%

LABELED COMPOUNDS

13C-2,3,7,8-TCDF 0.84 0.82 0.94 0.87 0.88 0.87 5.47 +/-35%

Analyst: \_\_\_\_\_ f Date: 3/4/09

6DFB - FORM VI-HR CDD-2  
TCDF INITIAL CALIBRATION ION ABUNDANCE RATIO SUMMARY  
HIGH RESOLUTION

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALC CASE NO.: EPA05-2008-01 TO NO.: SDG NO.: SDG E3YC9

GC COLUMN: DB225 ID: 0.25 (mm) INSTRUMENT ID: FAL1

INIT. CALIB. DATE(S): 13-JAN-09 13-JAN-09 13-JAN-09 13-JAN-09 13-JAN-09

INIT. CALIB. TIMES: 11:41:31 12:16:35 11:06:28 12:51:39 13:26:43

TARGET ANALYTES	SELECTED IONS	ION ABUNDANCE RATIO					FLAG	ION RATIO QC LIMITS
		CS1	CS2	CS3	CS4	CS5		
2,3,7,8-TCDF	304/306	0.78	0.75	0.78	0.75	0.75		0.65-0.89

LABELED COMPOUNDS								
13C-2,3,7,8-TCDF	316/318	0.80	0.80	0.79	0.79	0.79		0.65-0.89

INTERNAL STANDARDS								
13C-1,2,3,4-TCDF	316/318	0.79	0.79	0.80	0.79	0.79		0.65-0.89

Quality Control (QC) limits represent +/-15% window around the theoretical ion abundance ratio. The laboratory must flag any analyte in any calibration solution which does not meet the ion abundance ratio QC limit by placing an asterik in the flag column

Analyst:

Date: 3/11/09

5DFC - FORM V-HR CDD-3  
TCDF ANALYTICAL SEQUENCE SUMMARY  
HIGH RESOLUTION

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA08-2008-01 TO NO.: SDG NO.: SDG E3YC9

GC Column: DB225 ID: 0.25 (mm) INSTRUMENT ID: FAL1

INIT. CALIB. DATE (S) : 13-JAN-09

INIT. CALIB. TIMES: 11:41:31, 12:16:35, 11:06:28, 12:51:39, 13:26:43

THE ANALYTICAL SEQUENCE OF STANDARDS, SAMPLES, BLANKS, AND LABORATORY CONTROL SAMPLES (LCSS) IS AS FOLLOWS:

EPA SAMPLE NO.	LAB SAMPLE ID	LAB FILE ID	DATE ANALYZED	TIME ANALYZED
1613 CS3 080827J	ST011309A3	13JAN09A	13-JAN-09	11:06:28
1613 CS1 080827H	ST011309A1	13JAN09A	13-JAN-09	11:41:31
1613 CS2 080827I	ST011309A2	13JAN09A	13-JAN-09	12:16:35
1613 CS4 080827K	ST011309A4	13JAN09A	13-JAN-09	12:51:39
1613 CS5 080827L	ST011309A5	13JAN09A	13-JAN-09	13:26:43
ZZZZ	ZZZZ	13JAN09A	13-JAN-09	14:01:47

ANALYST: f DATE: 3/4/09

## Frontier Analytical Laboratory

Data Filename: 13JAN09A

Analyte: TCDFAL1-1-13-09

Cal: TCDFAL1-1-13-09

Name	RRF	S. D.	%RSD	S2 RRF#1	S3 RRF#2	S1 RRF#3	S4 RRF#4	S5 RRF#5	S5 RRF#6
2,3,7,8-TCDF	1.03	0.0651	6.35 %	1.12	0.97	1.06	1.00	0.98	0.98
13C-2,3,7,8-TCDF	0.87	0.0477	5.47 %	0.84	0.82	0.94	0.87	0.88	0.88
13C-1,2,3,4-TCDF	-	-	- %	-	-	-	-	-	-

Analyst: SLDate: 1/13/09

Run #1      Filename 13JAN09A      S: 2      Acquired: 13-JAN-09 11:41:31 Cal: TCDFFAL1-1-13-09  
Client ID: ST011309A1      Analyte:      FAL ID: 1613 CS1 080827H

Typ	Name	Amount	Resp	RA	RT	RF	RRF	
1	Unk	2,3,7,8-TCDF	0.50	9.58e+05	0.78 y	19:39	-	1.12 y
2	IS	13C-2,3,7,8-TCDF	100.00	1.71e+08	0.80 y	19:39	-	0.840 y
3	RS	13C-1,2,3,4-TCDF	100.00	2.04e+08	0.79 y	17:05	2.04e+06	- n

Analyst: E

Date: 1/13/09

Run #2      Filename 13JAN09A      S: 3      Acquired: 13-JAN-09 12:16:35      Cal: TCDFAL1-1-13-09  
Client ID: ST011309A2      Analyte:      FAL ID: 1613 CS2 080827I

Typ	Name	Amount	Resp	RA	RT	RF	RRF	
1	Unk	2,3,7,8-TCDF	2.00	2.82e+06	0.75 y	19:38	-	0.967 y
2	IS	13C-2,3,7,8-TCDF	100.00	1.46e+08	0.80 y	19:37	-	0.819 y
3	RS	13C-1,2,3,4-TCDF	100.00	1.78e+08	0.79 y	17:03	1.78e+06	- n

Analyst: E Date: 1/13/09

Run #3      Filename 13JAN09A      S: 1      Acquired: 13-JAN-09 11:06:28    Cal: TCDFAL1-1-13-09  
Client ID: ST011309A3      Analyte:      FAL ID: 1613 CS3 080827J

	Typ	Name	Amount	Resp	RA	RT	RF	RRF		
1	Unk	2,3,7,8-TCDF	10.00	1.70e+07	0.78	y	19:42	-	1.06	y
2	IS	13C-2,3,7,8-TCDF	100.00	1.60e+08	0.79	y	19:40	-	0.944	y
3	RS	13C-1,2,3,4-TCDF	100.00	1.69e+08	0.80	y	17:06	1.69e+06	-	n

Analyst: SC

Date: 4/13/09

Run #4      Filename 13JAN09A      S: 4      Acquired: 13-JAN-09 12:51:39    Cal: TCDFAL1-1-13-09  
Client ID: ST011309A4      Analyte:      FAL ID: 1613 CS4 080827K

	Typ	Name	Amount	Resp	RA	RT	RF	RRF		
1	Unk	2,3,7,8-TCDF	40.00	6.46e+07	0.75	y	19:39	-	1.00	y
2	IS	13C-2,3,7,8-TCDF	100.00	1.61e+08	0.79	y	19:37	-	0.872	y
3	RS	13C-1,2,3,4-TCDF	100.00	1.85e+08	0.79	y	17:04	1.85e+06	-	n

Analyst: E

Date: 1/13/09


Run #5      Filename 13JAN09A      S: 5      Acquired: 13-JAN-09 13:26:43 Cal: TCDFFAL1-1-13-09  
Client ID: ST011309A5      Analyte: TCDFFAL1-1-13-09      FAL ID: 1613 CS5 080827L

	Typ	Name	Amount	Resp	RA	RT	RF	RRF		
1	Unk	2,3,7,8-TCDF	200.00	3.34e+08	0.75	y	19:39	-	0.976	y
2	IS	13C-2,3,7,8-TCDF	100.00	1.71e+08	0.79	y	19:38	-	0.881	y
3	RS	13C-1,2,3,4-TCDF	100.00	1.94e+08	0.79	y	17:06	1.94e+06	-	n


Analyst: C

Date: 1/13/09

Peak Locate Examination:13-JAN-2009:11:05 File:13JAN09A  
Experiment:TCDF Function:1 Reference:PFK



Peak Locate Examination:13-JAN-2009:14:38 File:13JAN09A.RES\_CHECK  
Experiment:TCDF Function:1 Reference:PFK



File:13JAN09A Acq:13-JAN-2009 11:06:28 FAL-1 Frontier Analytical Laboratory  
303.9016 Expt:TCDF  
FAL ID:ST011309A3

100 %

1.8E6

1.6E6

1.4E6

1.2E6

1.0E6

8.9E5


7.1E5

5.3E5

3.6E5

1.8E5

0.0E0



File:13JAN09A Acq:13-JAN-2009 11:06:28 FAL-1 Frontier Analytical Laboratory  
317.9389 Expt:TCDF  
FAL ID:ST011309A3

100 %

1.6E7

1.5E7

1.3E7

1.1E7

9.9E6


8.2E6

6.6E6


4.9E6

3.3E6


1.6E6




File:13JAN09A Acq:13-JAN-2009 11:41:31 FAL-1 Frontier Analytical Laboratory  
 303,9016 S:2 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100.0,0.00%,F,F) Exp:TCDF  
 FAL ID:ST011309A1




File:13JAN09A Acq:13-JAN-2009 11:41:31 FAL-1 Frontier Analytical Laboratory  
 305,8987 S:2 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100.0,0.00%,F,F) Exp:TCDF  
 FAL ID:ST011309A1



File:13JAN09A Acq:13-JAN-2009 11:41:31 FAL-1 Frontier Analytical Laboratory  
 315,9419 S:2 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100.0,0.00%,F,F) Exp:TCDF  
 FAL ID:ST011309A1



File:13JAN09A Acq:13-JAN-2009 11:41:31 FAL-1 Frontier Analytical Laboratory  
 317,9389 S:2 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100.0,0.00%,F,F) Exp:TCDF  
 FAL ID:ST011309A1



File:13JAN09A Acq:13-JAN-2009 12:16:35 FAL-1 Frontier Analytical Laboratory

303.9016 S3 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100.0,0.00%,F,F) Exp:TCDF

FAL ID:ST011309A2

A1.21E6

2.6E5

1.3E5

0.0E0



File:13JAN09A Acq:13-JAN-2009 12:16:35 FAL-1 Frontier Analytical Laboratory

317.9389 S3 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100.0,0.00%,F,F) Exp:TCDF

FAL ID:ST011309A2

A7.86E7

2.0E4

1.1E4

0.0E0

H2.13E7

2.0E4

1.1E4

0.0E0

H1.60E7

2.0E4

1.1E4

0.0E0

A9.93E7

2.0E4

1.1E4

0.0E0

2.1E7

2.0E4

1.1E4

0.0E0

1.1E7

2.0E4

1.1E4

0.0E0

0.0E0

2.0E4

1.1E4

0.0E0

1.1E7

2.0E4

1.1E4

0.0E0

0.0E0

2.0E4

1.1E4

0.0E0

1.1E7

2.0E4

1.1E4

0.0E0

0.0E0

2.0E4

1.1E4

0.0E0

File:13JAN09A Acq:13-JAN-2009 11:06:28 FAL-1 Frontier Analytical Laboratory  
303.9016 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:TCDF  
FAL ID:ST011309A3

A1.78E7


A1.00E7

A1.48E7

4.1E6

2.1E6

0.0E0



File:13JAN09A Acq:13-JAN-2009 11:06:28 FAL-1 Frontier Analytical Laboratory  
315.9389 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:TCDF  
FAL ID:ST011309A3

A2.28E7

A1.29E7

A1.84E7

5.3E6

2.6E6

0.0E0

File:13JAN09A Acq:13-JAN-2009 11:06:28 FAL-1 Frontier Analytical Laboratory  
315.9419 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:TCDF  
FAL ID:ST011309A3

A7.52E7

A7.05E7

1.6E7

7.8E6

0.0E0

File:13JAN09A Acq:13-JAN-2009 11:06:28 FAL-1 Frontier Analytical Laboratory  
375.8364 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:TCDF  
FAL ID:ST011309A3

H1.92E7

H1.64E7

A8.92E7

1.9E7

9.6E6

0.0E0

File:13JAN09A Acq:13-JAN-2009 11:06:28 FAL-1 Frontier Analytical Laboratory  
375.8364 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:TCDF  
FAL ID:ST011309A3

H1.92E7

H1.64E7

A8.92E7

2.5E4

1.9E6

0.0E0

File:13JAN09A Acq:13-JAN-2009 11:06:28 FAL-1 Frontier Analytical Laboratory  
375.8364 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:TCDF  
FAL ID:ST011309A3

H1.92E7

H1.64E7

A8.92E7

1.9E7

9.6E6

0.0E0

File:13JAN09A Acq:13-JAN-2009 11:06:28 FAL-1 Frontier Analytical Laboratory  
375.8364 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:TCDF  
FAL ID:ST011309A3

H1.92E7

H1.64E7

A8.92E7

2.5E4

1.9E6

0.0E0

File:13JAN09A Acq:13-JAN-2009 12:51:39 FAL-1 Frontier Analytical Laboratory  
 303,9016 S:4 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100,0,0.00%,F,F) Exp:TCDF  
 FAL ID:ST011309A4




5.6E6  
2.8E6




7.5E6  
3.7E6



0.0E0  
0.0E0



1.8E7  
-8.9E6




2.3E7  
-1.1E7



0.0E0  
-0.0E0



2.1E4  
-1.1E4



-2.1E4  
-1.1E4

File:13JAN09A Acq:13-JAN-2009 13:26:43 FAL-1 Frontier Analytical Laboratory

303.9016 S:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:TCDF

FAL ID:ST011309A5

100 %

A1.43E8

3.0E7

-1.5E7

0.0E0

50 %

3.9E7

2.0E7

0.0E0

0 %

1.5E7

0.0E0

13:00 14:00 15:00 16:00 17:00 18:00 19:00 20:00 21:00 22:00 23:00 24:00 Time

File:13JAN09A Acq:13-JAN-2009 13:26:43 FAL-1 Frontier Analytical Laboratory

315.9419 S:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:TCDF

FAL ID:ST011309A5

100 %

A1.91E8

3.9E7

2.0E7

0.0E0

50 %

3.9E7

2.0E7

0.0E0

0 %

1.5E7

0.0E0

13:00 14:00 15:00 16:00 17:00 18:00 19:00 20:00 21:00 22:00 23:00 24:00 Time

File:13JAN09A Acq:13-JAN-2009 13:26:43 FAL-1 Frontier Analytical Laboratory

317.9389 S:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:TCDF

FAL ID:ST011309A5

100 %

A8.59E7

2.0E7

9.8E6

0.0E0

50 %

2.0E7

1.2E7

0.0E0

0 %

1.2E7

0.0E0

13:00 14:00 15:00 16:00 17:00 18:00 19:00 20:00 21:00 22:00 23:00 24:00 Time

File:13JAN09A Acq:13-JAN-2009 13:26:43 FAL-1 Frontier Analytical Laboratory

375.8364 S:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:TCDF

FAL ID:ST011309A5

100 %

H2.44E7

2.4E7

1.2E7

0.0E0

50 %

1.2E7

0.0E0

0 %

1.2E7

0.0E0

13:00 14:00 15:00 16:00 17:00 18:00 19:00 20:00 21:00 22:00 23:00 24:00 Time

File:13JAN09A Acq:13-JAN-2009 13:26:43 FAL-1 Frontier Analytical Laboratory

375.8364 S:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:TCDF

FAL ID:ST011309A5

100 %

A9.62E4

1.9E4

9.7E3

0.0E0

0 %

1.9E4

9.7E3

0.0E0

7DFA - Form VII-HR CDD-1  
CDD/CDF CONTINUING CALIBRATION SUMMARY  
HIGH RESOLUTION

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: SDG NO.: SDG E3YC9

GC COLUMN: DB5 ID: 0.25 (mm) INSTRUMENT ID: FAL3

LAB FILE ID: 24FEB09M Sam: 11 DATE ANALYZED: 23-FEB-09 TIME ANALYZED: 22:50:13

INIT. CALIB. TIMES: 16:15:11 INIT. CALIB. DATE(S): 5-FEB-09

TARGET ANALYTES	SELECTED IONS	MEAN				ION RATIO	ION FLAG	QC LIMITS
		RR/	RR/	%D	%D			
2,3,7,8-TCDD	320/322	1.07	1.04	2.69		0.79		0.65-0.89
2,3,7,8-TCDF	304/306	1.19	1.34	-11.4		0.72		0.65-0.89
1,2,3,7,8-PeCDF	340/342	1.02	1.03	-0.790		1.61		1.32-1.78
1,2,3,7,8-PeCDD	356/358	0.966	0.901	7.16		1.57		1.32-1.78
2,3,4,7,8-PeCDF	340/342	0.779	0.808	-3.64		1.60		1.32-1.78
1,2,3,4,7,8-HxCDF	374/376	1.44	1.35	6.20		1.28		1.05-1.43
1,2,3,6,7,8-HxCDF	374/376	1.32	1.25	5.72		1.28		1.05-1.43
1,2,3,4,7,8-HxCDD	390/392	1.37	1.43	-3.77		1.30		1.05-1.43
1,2,3,6,7,8-HxCDD	390/392	0.975	1.02	-4.34		1.25		1.05-1.43
1,2,3,7,8,9-HxCDD	390/392	1.08	1.05	2.36		1.25		1.05-1.43
2,3,4,6,7,8-HxCDF	374/376	1.58	1.53	3.52		1.27		1.05-1.43
1,2,3,7,8,9-HxCDF	374/376	1.22	1.19	2.60		1.27		1.05-1.43
1,2,3,4,6,7,8-HpCDF	408/410	1.17	1.10	6.07		1.03		0.88-1.20
1,2,3,4,6,7,8-HpCDD	424/426	0.884	0.987	-10.4		1.05		0.88-1.20
1,2,3,4,7,8,9-HpCDF	408/410	1.08	1.02	5.52		1.02		0.88-1.20
OCDD	458/460	0.884	0.940	-5.96		0.90		0.76-1.02
OCDF	442/444	0.747	0.763	-2.13		0.91		0.76-1.02

LABELED COMPOUNDS							
13C-2,3,7,8-TCDD	332/334	0.952	0.937	1.67		0.77	0.65-0.89
13C-1,2,3,7,8-PeCDD	368/370	0.902	0.754	19.6		1.78	1.32-1.78
13C-1,2,3,4,7,8-HxCDD	402/404	0.980	1.03	-4.66		1.30	1.05-1.43
13C-1,2,3,6,7,8-HxCDD	402/404	1.01	1.08	-6.05		1.29	1.05-1.43
13C-1,2,3,4,6,7,8-HpCDD	436/438	0.913	0.894	2.23		1.08	0.88-1.20
13C-OCDD	470/472	0.614	0.635	-3.30		0.93	0.76-1.02
13C-2,3,7,8-TCDF	316/318	0.889	0.863	3.10		0.85	0.65-0.89
13C-1,2,3,7,8-PeCDF	352/354	0.814	0.758	7.35		1.66	1.32-1.78
13C-2,3,4,7,8-PeCDF	352/354	0.844	0.776	8.73		1.67	1.32-1.78
13C-1,2,3,4,7,8-HxCDF	384/386	1.24	1.34	-7.60		0.53	0.43-0.59
13C-1,2,3,6,7,8-HxCDF	384/386	1.28	1.40	-8.23		0.53	0.43-0.59
13C-1,2,3,7,8,9-HxCDF	384/386	0.993	1.01	-1.29		0.54	0.43-0.59
13C-2,3,4,6,7,8-HxCDF	384/386	1.22	1.29	-5.62		0.54	0.43-0.59
13C-1,2,3,4,6,7,8-HpCDF	418/420	1.03	1.16	-11.2		0.44	0.37-0.51
13C-1,2,3,4,7,8,9-HpCDF	418/420	0.898	0.957	-6.12		0.44	0.37-0.51
13C-OCDF	454/456	0.843	0.947	-11.1		0.91	0.76-1.02

CLEAN-UP							
37Cl-2,3,7,8-TCDD	328/NA	0.847	0.879	-3.73		NA	NA

INTERNAL STANDARDS							
13C-1,2,3,4-TCDD	332/334	NA	NA	NA	NA	0.79	0.65-0.89
13C-1,2,3,4-TCDF	316/318	NA	NA	NA	NA	0.86	0.65-0.89
13C-1,2,3,7,8,9-HxCDD	402/404	NA	NA	NA	NA	1.29	1.05-1.43

The laboratory must flag any analyte which does not meet criteria

for Percent Difference (%D) or ion abundance ratio

by placing an asterisk in the appropriate flag column

Analyst: SC

Date: 3/3/09

7DFA - Form VII-HR CDD-1  
CDD/CDF CONTINUING CALIBRATION SUMMARY  
HIGH RESOLUTION

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: SDG NO.: SDG E3YC9

GC COLUMN: DB5 ID: 0.25 (mm) INSTRUMENT ID: FAL3

LAB FILE ID: 24FEB09M Sam: 18 DATE ANALYZED: 24-FEB-09 TIME ANALYZED: 05:17:02

INIT. CALIB. TIMES: 16:15:11 INIT. CALIB. DATE(S): 5-FEB-09

	SELECTED	MEAN			ION		
TARGET ANALYTES	IONS	RR/	RRF	%D	FLAG	RATIO	ION RATIO
2,3,7,8-TCDD	320/322	1.02	1.04	-2.17		0.81	0.65-0.89
2,3,7,8-TCDF	304/306	1.18	1.34	-12.1		0.71	0.65-0.89
1,2,3,7,8-PeCDF	340/342	1.02	1.03	-0.664		1.62	1.32-1.78
1,2,3,7,8-PeCDD	356/358	0.925	0.901	2.66		1.56	1.32-1.78
2,3,4,7,8-PeCDF	340/342	0.770	0.808	-4.71		1.60	1.32-1.78
1,2,3,4,7,8-HxCDF	374/376	1.44	1.35	6.20		1.26	1.05-1.43
1,2,3,6,7,8-HxCDF	374/376	1.31	1.25	4.86		1.25	1.05-1.43
1,2,3,4,7,8-HxCDD	390/392	1.39	1.43	-2.78		1.27	1.05-1.43
1,2,3,6,7,8-HxCDD	390/392	0.974	1.02	-4.44		1.29	1.05-1.43
1,2,3,7,8,9-HxCDD	390/392	1.10	1.05	4.37		1.29	1.05-1.43
2,3,4,6,7,8-HxCDF	374/376	1.56	1.53	2.02		1.25	1.05-1.43
1,2,3,7,8,9-HxCDF	374/376	1.22	1.19	3.12		1.28	1.05-1.43
1,2,3,4,6,7,8-HpCDF	408/410	1.16	1.10	5.50		1.02	0.88-1.20
1,2,3,4,6,7,8-HpCDD	424/426	0.892	0.987	-9.60		1.04	0.88-1.20
1,2,3,4,7,8,9-HpCDF	408/410	1.11	1.02	8.01		1.03	0.88-1.20
OCDD	458/460	0.878	0.940	-6.54		0.88	0.76-1.02
OCDF	442/444	0.746	0.763	-2.24		0.91	0.76-1.02

LABELED COMPOUNDS

13C-2,3,7,8-TCDD	332/334	0.963	0.937	2.81		0.78	0.65-0.89
13C-1,2,3,7,8-PeCDD	368/370	0.903	0.754	19.8		1.77	1.32-1.78
13C-1,2,3,4,7,8-HxCDD	402/404	0.953	1.03	-7.31		1.29	1.05-1.43
13C-1,2,3,6,7,8-HxCDD	402/404	1.02	1.08	-5.08		1.28	1.05-1.43
13C-1,2,3,4,6,7,8-HpCDD	436/438	0.858	0.894	-3.92		1.06	0.88-1.20
13C-OCDD	470/472	0.592	0.635	-6.80		0.93	0.76-1.02
13C-2,3,7,8-TCDF	316/318	0.890	0.863	3.14		0.85	0.65-0.89
13C-1,2,3,7,8-PeCDF	352/354	0.791	0.758	4.36		1.65	1.32-1.78
13C-2,3,4,7,8-PeCDF	352/354	0.810	0.776	4.45		1.66	1.32-1.78
13C-1,2,3,4,7,8-HxCDF	384/386	1.25	1.34	-6.95		0.53	0.43-0.59
13C-1,2,3,6,7,8-HxCDF	384/386	1.29	1.40	-7.53		0.54	0.43-0.59
13C-1,2,3,7,8,9-HxCDF	384/386	0.996	1.01	-0.990		0.54	0.43-0.59
13C-2,3,4,6,7,8-HxCDF	384/386	1.20	1.29	-6.87		0.53	0.43-0.59
13C-1,2,3,4,6,7,8-HpCDF	418/420	1.02	1.16	-11.6		0.44	0.37-0.51
13C-1,2,3,4,7,8,9-HpCDF	418/420	0.853	0.957	-10.9		0.45	0.37-0.51
13C-OCDF	454/456	0.802	0.947	-15.4		0.92	0.76-1.02

CLEAN-UP

37Cl-2,3,7,8-TCDD	328/NA	0.877	0.879	-0.271	NA	NA	NA
-------------------	--------	-------	-------	--------	----	----	----

INTERNAL STANDARDS

13C-1,2,3,4-TCDD	332/334	NA	NA	NA	NA	0.79	0.65-0.89
13C-1,2,3,4-TCDF	316/318	NA	NA	NA	NA	0.84	0.65-0.89
13C-1,2,3,7,8,9-HxCDD	402/404	NA	NA	NA	NA	1.28	1.05-1.43

The laboratory must flag any analyte which does not meet criteria for Percent Difference (%D) or ion abundance ratio by placing an asterik in the appropriate flag column

Analyst: AC

Date: 3/6/09

7DFA - Form VII-HR CDD-1  
 TCDF CONTINUING CALIBRATION SUMMARY  
 HIGH RESOLUTION

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: SDG NO.: SDG E3YC9

GC COLUMN: DB225 ID: 0.25 (mm) INSTRUMENT ID: FAL1

LAB FILE ID: 04MAR09A DATE ANALYZED: 4-MAR-09 TIME ANALYZED: 09:46:39

INIT. CALIB. TIMES: 11:41:31 INIT. CALIB. DATE(S): 13-JAN-09

TARGET ANALYTES	SELECTED IONS	MEAN			ION		
		RR/ RRF	RR/ RRF	%D	FLAG	ION RATIO	QC LIMITS
2,3,7,8-TCDF	304/306	0.918	1.03	-10.6	0.78	0.65-0.89	

LABELED COMPOUNDS							
13C-2,3,7,8-TCDF	316/318	0.909	0.871	4.35	0.79	0.65-0.89	

CLEAN-UP							
	328/NA	NA	NA	NA	NA	NA	NA

INTERNAL STANDARDS							
13C-1,2,3,4-TCDF	316/318	NA	NA	NA	NA	0.79	0.65-0.89

The laboratory must flag any analyte which does not meet criteria  
 for Percent Difference (%D) or ion abundance ratio  
 by placing an asterik in the appropriate flag column

Analyst: J

Date: 3/4/09

7DFA - Form VII-HR CDD-1  
 TCDF CONTINUING CALIBRATION SUMMARY  
 HIGH RESOLUTION

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: SDG NO.: SDG E3YC9

GC COLUMN: DB225 ID: 0.25 (mm) INSTRUMENT ID: FAL1

LAB FILE ID: 04MAR09A DATE ANALYZED: 4-MAR-09 TIME ANALYZED: 12:09:12

INIT. CALIB. TIMES: 11:41:31 INIT. CALIB. DATE(S): 13-JAN-09

TARGET ANALYTES	SELECTED IONS	MEAN				ION		
		RR/	RRF	%D	%D	ION RATIO	FLAG	QC LIMITS
2,3,7,8-TCDF	304/306	0.841	1.03	-18.0	0.79		0.65-0.89	

LABELED COMPOUNDS	316/318	0.811	0.871	-6.93	0.78	0.65-0.89
13C-2,3,7,8-TCDF						

CLEAN-UP	328/NA	NA	NA	NA	NA	NA	NA
INTERNAL STANDARDS							
13C-1,2,3,4-TCDF	316/318	NA	NA	NA	NA	0.78	0.65-0.89

The laboratory must flag any analyte which does not meet criteria  
 for Percent Difference (%D) or ion abundance ratio  
 by placing an asterik in the appropriate flag column

Analyst: 

Date: 3/4/09

7DFB - Form VII-HR CDD-2  
 CDD/CDF CONTINUING CALIBRATION RETENTION TIME SUMMARY  
 HIGH RESOLUTION

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: SDG NO.: SDG E3YC9

GC COLUMN: DB5 ID: 0.25 (mm) INSTRUMENT ID: FAL3

DATE ANALYZED: 23-FEB-09 TIME ANALYZED: 22:50:13

LAB FILE ID: 24FEB09M Sam: 11

INIT. CALIB. TIMES: 16:15:11 INIT. CALIB. DATE(S): 5-FEB-09

TARGET ANALYTES	RRT	RT
2,3,7,8-TCDD	1.001	27:32
2,3,7,8-TCDF	1.001	26:46
1,2,3,7,8-PeCDF	1.001	31:37
1,2,3,7,8-PeCDD	1.001	33:20
2,3,4,7,8-PeCDF	1.000	32:55
1,2,3,4,7,8-HxCDF	1.001	37:18
1,2,3,6,7,8-HxCDF	1.001	37:31
1,2,3,4,7,8-HxCDD	1.001	38:42
1,2,3,6,7,8-HxCDD	1.001	38:52
1,2,3,7,8,9-HxCDD	1.012	39:18
2,3,4,6,7,8-HxCDF	1.000	38:26
1,2,3,7,8,9-HxCDF	1.000	39:52
1,2,3,4,6,7,8-HpCDF	1.001	42:24
1,2,3,4,6,7,8-HpCDD	1.000	44:17
1,2,3,4,7,8,9-HpCDF	1.001	45:13
OCDD	1.001	49:52
OCDF	1.001	50:15

LABELED COMPOUNDS

13C-2,3,7,8-TCDD	1.022	27:30
13C-1,2,3,7,8-PeCDD	1.238	33:19
13C-1,2,3,4,7,8-HxCDD	0.985	38:40
13C-1,2,3,6,7,8-HxCDD	0.989	38:50
13C-1,2,3,4,6,7,8-HpCDD	1.127	44:17
13C-OCDD	1.269	49:51
13C-2,3,7,8-TCDF	0.994	26:45
13C-1,2,3,7,8-PeCDF	1.173	31:35
13C-2,3,4,7,8-PeCDF	1.223	32:55
13C-1,2,3,4,7,8-HxCDF	0.949	37:17
13C-1,2,3,6,7,8-HxCDF	0.954	37:29
13C-1,2,3,7,8,9-HxCDF	1.015	39:52
13C-2,3,4,6,7,8-HxCDF	0.978	38:26
13C-1,2,3,4,6,7,8-HpCDF	1.079	42:22
13C-1,2,3,4,7,8,9-HpCDF	1.151	45:12
13C-OCDF	1.279	50:14

CLEAN-UP STANDARD

37Cl-2,3,7,8-TCDD	1.023	27:32
-------------------	-------	-------

INTERNAL STANDARD

13C-1,2,3,4-TCDD	NA	26:55
13C-1,2,3,4-TCDF	NA	25:40
13C-1,2,3,7,8,9-HxCDD	NA	39:17

Analyst: C

Date: 3/2/09

7DFB - Form VII-HR CDD-2  
 CDD/CDF CONTINUING CALIBRATION RETENTION TIME SUMMARY  
 HIGH RESOLUTION

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: SDG NO.: SDG E3YC9

GC COLUMN: DB5 ID: 0.25 (mm) INSTRUMENT ID: FAL3

DATE ANALYZED: 24-FEB-09 TIME ANALYZED: 05:17:02

LAB FILE ID: 24FEB09M Sam: 18

INIT. CALIB. TIMES: 16:15:11 INIT. CALIB. DATE(S): 5-FEB-09

TARGET ANALYTES	RRT	RT
2,3,7,8-TCDD	1.001	27:31
2,3,7,8-TCDF	1.001	26:46
1,2,3,7,8-PeCDF	1.001	31:36
1,2,3,7,8-PeCDD	1.001	33:20
2,3,4,7,8-PeCDF	1.001	32:55
1,2,3,4,7,8-HxCDF	1.001	37:18
1,2,3,6,7,8-HxCDF	1.000	37:29
1,2,3,4,7,8-HxCDD	1.001	38:41
1,2,3,6,7,8-HxCDD	1.000	38:51
1,2,3,7,8,9-HxCDD	1.012	39:18
2,3,4,6,7,8-HxCDF	1.001	38:26
1,2,3,7,8,9-HxCDF	1.001	39:53
1,2,3,4,6,7,8-HpCDF	1.001	42:22
1,2,3,4,6,7,8-HpCDD	1.001	44:17
1,2,3,4,7,8,9-HpCDF	1.001	45:13
OCDD	1.001	49:52
OCDF	1.001	50:15

LABELED COMPOUNDS

13C-2,3,7,8-TCDD	1.022	27:30
13C-1,2,3,7,8-PeCDD	1.238	33:18
13C-1,2,3,4,7,8-HxCDD	0.985	38:40
13C-1,2,3,6,7,8-HxCDD	0.989	38:50
13C-1,2,3,4,6,7,8-HpCDD	1.127	44:15
13C-OCDD	1.269	49:50
13C-2,3,7,8-TCDF	0.994	26:45
13C-1,2,3,7,8-PeCDF	1.174	31:34
13C-2,3,4,7,8-PeCDF	1.222	32:53
13C-1,2,3,4,7,8-HxCDF	0.949	37:16
13C-1,2,3,6,7,8-HxCDF	0.954	37:29
13C-1,2,3,7,8,9-HxCDF	1.015	39:51
13C-2,3,4,6,7,8-HxCDF	0.978	38:24
13C-1,2,3,4,6,7,8-HpCDF	1.078	42:21
13C-1,2,3,4,7,8,9-HpCDF	1.151	45:11
13C-OCDF	1.279	50:13

CLEAN-UP STANDARD

37Cl-2,3,7,8-TCDD	1.023	27:31
-------------------	-------	-------

INTERNAL STANDARD

13C-1,2,3,4-TCDD	NA	26:54
13C-1,2,3,4-TCDF	NA	25:39
13C-1,2,3,7,8,9-HxCDD	NA	39:16

Analyst: AC

Date: 3/9/09

7DFB - Form VII-HR CDD-2  
TCDF CONTINUING CALIBRATION RETENTION TIME SUMMARY  
HIGH RESOLUTION

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: SDG NO.: SDG E3YC9

GC COLUMN: DB225 ID: 0.25 (mm) INSTRUMENT ID: FAL1

DATE ANALYZED: 4-MAR-09 TIME ANALYZED: 09:46:39

LAB FILE ID: 04MAR09A

INIT. CALIB. TIMES: 11:41:31 INIT. CALIB. DATE(S): 13-JAN-09

TARGET ANALYTES	RRT	RT
2,3,7,8-TCDF	1.001	19:34

LABELED COMPOUNDS  
13C-2,3,7,8-TCDF NA 19:33

CLEAN-UP STANDARD

INTERNAL STANDARD  
13C-1,2,3,4-TCDF NA 16:59

Analyst: J

Date: 3/4/09

7DFB - Form VII-HR CDD-2  
TCDF CONTINUING CALIBRATION RETENTION TIME SUMMARY  
HIGH RESOLUTION

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: SDG NO.: SDG E3YC9

GC COLUMN: DB225 ID: 0.25 (mm) INSTRUMENT ID: FAL1

DATE ANALYZED: 4-MAR-09 TIME ANALYZED: 12:09:12

LAB FILE ID: 04MAR09A

INIT. CALIB. TIMES: 11:41:31 INIT. CALIB. DATE(S): 13-JAN-09

TARGET ANALYTES	RRT	RT
2,3,7,8-TCDF	1.002	19:29

LABELED COMPOUNDS  
13C-2,3,7,8-TCDF NA 19:27


CLEAN-UP STANDARD

INTERNAL STANDARD  
13C-1,2,3,4-TCDF NA 16:54


Analyst: J

Date: 3/4/09


Peak Locate Examination:23-FEB-2009:13:35 File:24FEB09M  
Experiment:PCDD Function:1 Reference:PFK




Peak Locate Examination:23-FEB-2009:13:35 File:24FEB09M  
Experiment:PCDD Function:2 Reference:PFK




Peak Locate Examination:23-FEB-2009:13:36 File:24FEB09M  
Experiment:PCDD Function:3 Reference:PFK




Peak Locate Examination:23-FEB-2009:13:36 File:24FEB09M  
Experiment:PCDD Function:4 Reference:PFK




Peak Locate Examination:23-FEB-2009:13:36 File:24FEB09M  
Experiment:PCDD Function:5 Reference:PFK




Peak Locate Examination:24-FEB-2009:07:47 File:24FEB09M\_RES\_CHECK  
Experiment:PCDD Function:1 Reference:PFK




Peak Locate Examination:24-FEB-2009:07:47 File:24FEB09M.RES\_CHECK  
Experiment:PCDD Function:2 Reference:PFK




Peak Locate Examination:24-FEB-2009:07:48 File:24FEB09M\_RES\_CHECK  
Experiment:PCDD Function:3 Reference:PFK



Peak Locate Examination:24-FEB-2009:07:48 File:24FEB09M\_RES\_CHECK  
Experiment:pCDD Function:4 Reference:PFK



Peak Locate Examination:24-FEB-2009:07:49 File:24FEB09M\_RES\_CHECK  
Experiment:PCDD Function:5 Reference:PFK



File:24FEB09M #1-390 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima

319.8965 S:11 Exp:PCDD

Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory

100 %

3.1E5

2.8E5

2.5E5

2.2E5

1.9E5

1.6E5

1.2E5


9.4E4

6.2E4

3.1E4

0.0E0

10.5/  
6.9 = 15.3%



100 %


3.4E5

1.0E6




6.8E5

3.4E5


0.0E0




File:24FEB09M #1-390 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
319.8965 S:11 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-425 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
355.8546 S:11 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-425 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
357.8517 S:11 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-425 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
367.8949 S:11 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-425 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
369.8919 S:11 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-425 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
366.9792 S:11 F:2 Exp:PCDD  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-464 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
389.8156 S:11 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-464 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
401.8559 S:11 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC







File:24FEB09M #1-464 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
403.8530 S:11 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-464 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
380.9760 S:11 F:3 Exp:PCDD  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-541 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
423.7767 S:11 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,R,F) Exp:PCDD  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-347 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
457.7377 S:11 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,FF) Exp:PCDD  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-347 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
459.7348 S:11 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,FF) Exp:PCDD  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC

A4.70E6  
49:52  
A5.23E6



File:24FEB09M #1-347 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
469.7780 S:11 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,FF) Exp:PCDD  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC

49.51  
A1.08E7




File:24FEB09M #1-347 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
471.7750 S:11 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,FF) Exp:PCDD  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC

49:51  
A1.16E7




File:24FEB09M #1-347 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
454.9728 S:11 F:5 Exp:PCDD  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC


100%  
48:13 48:35 48:58 49:22 49:53 50:20 50:43 50:59 51:17 51:46 52:03




File:24FEB09M #1-390 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
303.9016 S:11 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-390 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
305.8987 S:11 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC







File:24FEB09M #1-390 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
315.9419 S:11 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-390 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
317.9389 S:11 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-390 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
339.8597 S:11 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-425 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
339.8597 S:11 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-425 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
351.9000 S:11 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-425 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
353.8970 S:11 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-425 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
409.7974 S:11 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC

File:24FEB09M #1-464 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
373.8207 S:11 F:3 BSUB(10000,15,-3,0) PKD(5,5,3,0.10%,100.0,0.00%,F,F)  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-464 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
375.8178 S:11 F:3 BSUB(10000,15,-3,0) PKD(5,5,3,0.10%,100.0,0.00%,F,F)  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-464 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
383.8639 S:11 F:3 BSUB(10000,15,-3,0) PKD(5,5,3,0.10%,100.0,0.00%,F,F)  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-464 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
385.8610 S:11 F:3 BSUB(10000,15,-3,0) PKD(5,5,3,0.10%,100.0,0.00%,F,F)  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-464 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
445.7555 S:11 F:3 BSUB(10000,15,-3,0) PKD(5,5,3,0.10%,100.0,0.00%,F,F)  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-541 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
407.7818 S:11 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-541 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
417.8253 S:11 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-541 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
419.8220 S:11 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC







File:24FEB09M #1-541 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
419.7165 S:11 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-541 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
419.7165 S:11 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-347 Acq:23-FEB-2009 22:50:13 GC EI+ Voltage SIR Autospec-Ultima  
441.7428 S:11 F:5 BSUB(10000,15,-3,0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST024109M3 File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-390 Acq:24-FEB-2009 05:17:02 GC EI+ Voltage SIR Autospec-Ultima

319,8965 S:18 Exp:PCDD

Sample Text:ST022409M4 File Text:Frontier Analytical Laboratory

100 %

2.5E5

2.2E5

2.0E5

1.7E5

1.5E5

1.2E5

9.8E4

7.4E4

4.9E4

2.5E4


0.0E0

27:23

27:31

28:27

$$q_{15}/q_{69} = 13.8$$



File:24FEB09M #1-390 Acq:24-FEB-2009 05:17:02 GC EI+ Voltage SIR Autospec-Ultima

333,9339 S:18 Exp:PCDD

Sample Text:ST022409M4 File Text:Frontier Analytical Laboratory

100 %

2.7E6

2.4E6

2.1E6

1.9E6

1.6E6

1.3E6


1.1E6

8.1E5


5.4E5

2.7E5


0.0E0




File:24FEB09M #1-390 Acq:24-FEB-2009 05:17:02 GC EI+ Voltage SIR Autospec-Ultima  
319.8965 S:18 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST022409M4 File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-390 Acq:24-FEB-2009 05:17:02 GC EI+ Voltage SIR Autospec-Ultima  
321.8936 S:18 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST022409M4 File Text:Frontier Analytical Laboratory Analyst ID:TC






File:24FEB09M #1-390 Acq:24-FEB-2009 05:17:02 GC EI+ Voltage SIR Autospec-Ultima  
333.9339 S:18 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST022409M4 File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-390 Acq:24-FEB-2009 05:17:02 GC EI+ Voltage SIR Autospec-Ultima  
333.9339 S:18 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST022409M4 File Text:Frontier Analytical Laboratory Analyst ID:TC


File:24FEB09M #1-425 Acq:24-FEB-2009 05:17:02 GC EI+ Voltage SIR Autospec-Ultima  
355.8546 S:18 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST022409M4 File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-464 Acq:24-FEB-2009 05:17:02 GC EI+ Voltage SIR Autospec-Ultima  
389.8156 S:18 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST022409M4 File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-464 Acq:24-FEB-2009 05:17:02 GC EI+ Voltage SIR Autospec-Ultima  
391.8127 S:18 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST022409M4 File Text:Frontier Analytical Laboratory Analyst ID:TC






File:24FEB09M #1-464 Acq:24-FEB-2009 05:17:02 GC EI+ Voltage SIR Autospec-Ultima  
401.8559 S:18 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST022409M4 File Text:Frontier Analytical Laboratory Analyst ID:TC







File:24FEB09M #1-464 Acq:24-FEB-2009 05:17:02 GC EI+ Voltage SIR Autospec-Ultima  
380.9760 S:18 F:3 Exp:PCDD  
Sample Text:ST022409M4 File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-541 Acq:24-FEB-2009 05:17:02 GC EI+ Voltage SIR Autospec-Ultima  
423.7767 S:18 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST022409M4 File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-347 Acq:24-FEB-2009 05:17:02 GC EI+ Voltage SIR Autospec-Ultima  
457.7377 S:18 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST022409M4 File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-390 Acq:24-FEB-2009 05:17:02 GC EI+ Voltage SIR Autospec-Ultima  
303.9016 S:18 BSUB(10000,15,-3,0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST022409M4 File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-390 Acq:24-FEB-2009 05:17:02 GC EI+ Voltage SIR Autospec-Ultima  
305.8987 S:18 BSUB(10000,15,-3,0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST022409M4 File Text:Frontier Analytical Laboratory Analyst ID:TC






File:24FEB09M #1-390 Acq:24-FEB-2009 05:17:02 GC EI+ Voltage SIR Autospec-Ultima  
315.9419 S:18 BSUB(10000,15,-3,0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST022409M4 File Text:Frontier Analytical Laboratory Analyst ID:TC







File:24FEB09M #1-390 Acq:24-FEB-2009 05:17:02 GC EI+ Voltage SIR Autospec-Ultima  
375.8364 S:18 BSUB(10000,15,-3,0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST022409M4 File Text:Frontier Analytical Laboratory Analyst ID:TC






File:24FEB09M #1-390 Acq:24-FEB-2009 05:17:02 GC EI+ Voltage SIR Autospec-Ultima  
339.8597 S:18 BSUB(10000,15,-3,0) PKD(5,5,3,0,10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST022409M4 File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-425 Acq:24-FEB-2009 05:17:02 GC EI+ Voltage SIR Autospec-Ultima  
339.8597 S:18 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST022409M4 File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-464 Acq:24-FEB-2009 05:17:02 GC EI+ Voltage SIR Autospec-Ultima  
373.8207 S:18 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST022409M4 File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-541 Acq:24-FEB-2009 05:17:02 GC EI+ Voltage SIR Autospec-Ultima  
407.7818 S:18 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST022409M4 File Text:Frontier Analytical Laboratory Analyst ID:TC

100 %  
A4.13E6  
42:22

A3.30E6  
45:13

8.1E5  
4.1E5  
0.0E0




File:24FEB09M #1-541 Acq:24-FEB-2009 05:17:02 GC EI+ Voltage SIR Autospec-Ultima  
417.8253 S:18 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST022409M4 File Text:Frontier Analytical Laboratory Analyst ID:TC

100 %  
A4.31E6  
42:21

A3.21E6  
45:13

7.9E5  
3.9E5  
0.0E0




File:24FEB09M #1-541 Acq:24-FEB-2009 05:17:02 GC EI+ Voltage SIR Autospec-Ultima  
419.8220 S:18 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST022409M4 File Text:Frontier Analytical Laboratory Analyst ID:TC

100 %  
A4.31E6  
42:21

A3.65E6  
45:11

8.3E5  
4.1E5  
0.0E0



File:24FEB09M #1-541 Acq:24-FEB-2009 05:17:02 GC EI+ Voltage SIR Autospec-Ultima  
479.7165 S:18 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:ST022409M4 File Text:Frontier Analytical Laboratory Analyst ID:TC

100 %  
A1.14E3  
42:21

A1.04E3  
43:02

A986.42  
43:33

A1.05E3  
44:02

A1.61E3  
44:48

A1.26E3  
45:37

A818.62  
46:08

A1.11E3  
46:35

A1.23E3  
47:01

A351.53  
47:30

A2.49E3  
44:38

A1.14E3  
42:00

A1.04E3  
43:00

A986.42  
44:00

A1.05E3  
45:00

A1.61E3  
46:00

A1.26E3  
47:00

A818.62  
48:00

1.3E3  
6.6E2  
0.0E0

Time

File:24FEB09M #1-347 Acq:24-FEB-2009 05:17:02 GC EI+ Voltage SIR Autospec-Ultima  
441.7428 S:18 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST022409M4 File Text:Frontier Analytical Laboratory Analyst ID:TC

100 %  
50  
0  
48:00 49:00 50:00 51:00 52:00 Time  
50:15 A3.93E6  
6.0E5  
3.0E5  
0.0E0

File:24FEB09M #1-347 Acq:24-FEB-2009 05:17:02 GC EI+ Voltage SIR Autospec-Ultima  
443.7398 S:18 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST022409M4 File Text:Frontier Analytical Laboratory Analyst ID:TC

100 %  
50  
0  
48:00 49:00 50:00 51:00 52:00 Time  
50:15 A4.32E6  
6.5E5  
3.3E5  
0.0E0

File:24FEB09M #1-347 Acq:24-FEB-2009 05:17:02 GC EI+ Voltage SIR Autospec-Ultima  
453.7831 S:18 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST022409M4 File Text:Frontier Analytical Laboratory Analyst ID:TC

100 %  
50  
0  
48:00 49:00 50:00 51:00 52:00 Time  
50:13 A1.06E7  
1.6E6  
7.8E5  
0.0E0


File:24FEB09M #1-347 Acq:24-FEB-2009 05:17:02 GC EI+ Voltage SIR Autospec-Ultima  
455.7801 S:18 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST022409M4 File Text:Frontier Analytical Laboratory Analyst ID:TC

100 %  
50  
0  
48:00 49:00 50:00 51:00 52:00 Time  
50:13 A1.15E7  
1.7E6  
8.6E5  
0.0E0


File:24FEB09M #1-347 Acq:24-FEB-2009 05:17:02 GC EI+ Voltage SIR Autospec-Ultima  
513.6775 S:18 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:ST022409M4 File Text:Frontier Analytical Laboratory Analyst ID:TC

100 %  
50  
0  
48:00 49:00 50:00 51:00 52:00 Time  
50:31 A718.79  
52:13 A918.44  
6.1E2  
8.6E5  
0.0E0  
48:36 A471.02  
48:57 A673.35  
50:16 A390.08  
50:58 A578.99  
51:29 A369.88  
51:36 A559.24  
52:20 A4356.26  
52:13 A4356.26  
0.0E0

Peak Locate Examination: 4-MAR-2009:09:46 File:04MAR09A  
 Experiment:TCDF Function:1 Reference:PFK



Peak Locate Examination: 4-MAR-2009:12:50 File:04MAR09A\_RES\_CHECK  
 Experiment:TCDF Function:1 Reference:PFK



File:04MAR09A Acc: 4-MAR-2009 09:46:39 FAL-1 Frontier Analytical Laboratory

303,9016 Exp:TCDF

FAL ID:ST030409A1

100 %

1.7E6

1.5E6

1.4E6

1.2E6

1.0E6

8.5E5

6.8E5

5.1E5

3.4E5

1.7E5

0.0E0


$$S/A \approx 7.05\%$$

19:44

19:39

19:34

19:33



File:04MAR09A #1-1388 Acq: 4-MAR-2009 09:46:39 GC EI+ Voltage SIR 70S  
303.9016 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100.0,0.00%,F,F) Exp:TCDF  
Sample Text:ST030409A1 File Text:Frontier Analytical Laboratory Analyst ID:BS

13:37  
A1.61E7  
17:41  
A1.67E7  
22:07  
A1.34E7

17:40  
A2.12E7

18:48  
A1.74E6  
19:45  
A1.17E7

22:07  
A1.70E7  
5.0E6  
2.5E6  
0.0E0

1.7E7  
8.5E6  
0.0E0

2.2E7  
1.1E7  
0.0E0

File:04MAR09A #1-1388 Acq: 4-MAR-2009 09:46:39 GC EI+ Voltage SIR 70S  
315.9419 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100.0,0.00%,F,F) Exp:TCDF  
Sample Text:ST030409A1 File Text:Frontier Analytical Laboratory Analyst ID:BS

13:36  
A2.09E7  
17:40  
A2.12E7

18:48  
A1.74E6  
19:45  
A1.17E7

22:07  
A1.70E7  
5.0E6  
2.5E6  
0.0E0

1.7E7  
8.5E6  
0.0E0

File:04MAR09A #1-1388 Acq: 4-MAR-2009 09:46:39 GC EI+ Voltage SIR 70S  
317.9389 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100.0,0.00%,F,F) Exp:TCDF  
Sample Text:ST030409A1 File Text:Frontier Analytical Laboratory Analyst ID:BS

100 %  
50

16:59  
A7.86E7  
19:33  
A7.14E7

16:59  
A9.99E7  
19:33  
A9.08E7

20:14  
A2.60E4  
21:29  
A2.94E4  
22:21  
A5.64E4  
23:27  
A1.49E4  
1.4E4  
1.1E7  
0.0E0

File:04MAR09A #1-1388 Acq: 4-MAR-2009 09:46:39 GC EI+ Voltage SIR 70S  
375.8364 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100.0,0.00%,F,F) Exp:TCDF  
Sample Text:ST030409A1 File Text:Frontier Analytical Laboratory Analyst ID:BS


100 %  
50

13:49  
A2.80E4  
15:17  
A1.65E4  
16:42  
A3.00E4  
16:58  
A2.36E4





20:14  
A2.60E4  
21:29  
A2.94E4  
22:21  
A5.64E4  
23:27  
A1.49E4  
1.4E4  
1.1E7  
0.0E0

7.0E3  
0.0E0

File:04MAR09A Acq: 4-MAR-2009 12:09:12 FAL-1 Frontier Analytical Laboratory  
303.9016 S:5 Exp:TCDF  
FAL ID:ST030409A2  
100 %



File:04MAR09A #1-1388 Acq: 4-MAR-2009 12:09:12 GC EI+ Voltage SIR 70S  
303.9016 S:5 BSUB(10000,15,-3,0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:TCDF  
Sample Text:ST030409A2 File Text:Frontier Analytical Laboratory Analyst ID:BS



FAL ID: ST024109M3      Filename: 24FEB09M      Sam:11      Acquired: 23-FEB-09 22:50:13      ICAL: PCDDFAL3-2-5-09

Client ID: 1613 CS3 (080827J)      ConCal: ST022409M3      EndCal: ST022409M4

Results: 5340		GC Column: DB5		Amount: 1.000		NATO 1989 Tox:		102			
Name	Resp	RA	RT	RRF	Conc	Qual	Fac	Noise-1	Noise-2	DL	
2,3,7,8-TCDD	2.50e+06	0.79	y	27:32	1.04	10.3	2.50	-	-	*	
1,2,3,7,8-PeCDD	1.07e+07	1.57	y	33:20	0.90	53.6	2.50	-	-	*	
1,2,3,4,7,8-HxCDD	1.23e+07	1.30	y	38:42	1.43	48.1	2.50	-	-	*	
1,2,3,6,7,8-HxCDD	9.03e+06	1.25	y	38:52	1.02	47.8	2.50	-	-	*	
1,2,3,7,8,9-HxCDD	9.81e+06	1.25	y	39:18	1.05	51.2	2.50	-	-	*	
1,2,3,4,6,7,8-HpCDD	7.39e+06	1.05	y	44:17	0.99	44.8	2.50	-	-	*	
OCDD	9.93e+06	0.90	y	49:52	0.94	94.0	2.50	-	-	*	
2,3,7,8-TCDF	4.65e+06	0.72	y	26:46	1.34	8.86	2.50	-	-	*	
1,2,3,7,8-PeCDF	1.83e+07	1.61	y	31:37	1.03	49.6	2.50	-	-	*	
2,3,4,7,8-PeCDF	1.45e+07	1.60	y	32:55	0.81	48.2	2.50	-	-	*	
1,2,3,4,7,8-HxCDF	1.63e+07	1.28	y	37:18	1.35	53.1	2.50	-	-	*	
1,2,3,6,7,8-HxCDF	1.55e+07	1.28	y	37:31	1.25	52.9	2.50	-	-	*	
2,3,4,6,7,8-HxCDF	1.76e+07	1.27	y	38:26	1.53	51.8	2.50	-	-	*	
1,2,3,7,8,9-HxCDF	1.11e+07	1.27	y	39:52	1.19	51.3	2.50	-	-	*	
1,2,3,4,6,7,8-HpCDF	1.09e+07	1.03	y	42:24	1.10	53.0	2.50	-	-	*	
1,2,3,4,7,8,9-HpCDF	8.87e+06	1.02	y	45:13	1.02	52.8	2.50	-	-	*	
OCDF	1.15e+07	0.91	y	50:15	0.76	97.9	2.50	-	-	*	
										Rec	
13C-2,3,7,8-TCDD	2.34e+07	0.77	y	27:30	0.94	102				102	
13C-1,2,3,7,8-PeCDD	2.22e+07	1.78	y	33:19	0.75	120				120	
13C-1,2,3,4,7,8-HxCDD	1.79e+07	1.30	y	38:40	1.03	95.3				95.3	
13C-1,2,3,6,7,8-HxCDD	1.85e+07	1.29	y	38:50	1.08	94.0				94.0	
13C-1,2,3,4,6,7,8-HpCDD	1.67e+07	1.08	y	44:17	0.89	102				102	
13C-OCDD	2.25e+07	0.93	y	49:51	0.64	193				96.7	
13C-2,3,7,8-TCDF	3.92e+07	0.85	y	26:45	0.86	103				103	
13C-1,2,3,7,8-PeCDF	3.59e+07	1.66	y	31:35	0.76	107				107	
13C-2,3,4,7,8-PeCDF	3.72e+07	1.67	y	32:55	0.78	109				109	
13C-1,2,3,4,7,8-HxCDF	2.27e+07	0.53	y	37:17	1.34	92.4				92.4	
13C-1,2,3,6,7,8-HxCDF	2.35e+07	0.53	y	37:29	1.40	91.8				91.8	
13C-2,3,4,6,7,8-HxCDF	2.22e+07	0.54	y	38:26	1.29	94.4				94.4	
13C-1,2,3,7,8,9-HxCDF	1.82e+07	0.54	y	39:52	1.01	98.7				98.7	
13C-1,2,3,4,6,7,8-HpCDF	1.88e+07	0.44	y	42:22	1.16	88.8				88.8	
13C-1,2,3,4,7,8,9-HpCDF	1.64e+07	0.44	y	45:12	0.96	93.9				93.9	
13C-OCDF	3.08e+07	0.91	y	50:14	0.95	178				88.9	
37Cl-2,3,7,8-TCDD	2.08e+06			27:32	0.88	9.63				96.3	
13C-1,2,3,4-TCDD	2.46e+07	0.79	y	26:55	-	66.1					
13C-1,2,3,4-TCDF	4.41e+07	0.86	y	25:40	-	73.6					
13C-1,2,3,7,8,9-HxCDD	1.83e+07	1.29	y	39:17	-	91.0					
Total Tetra-Dioxins	1.28e+07			23:03	1.04	52.8	2.50	-	-	*	
Total Penta-Dioxins	3.50e+07			30:21	0.90	175	2.50	-	-	*	
Total Hexa-Dioxins	4.61e+07			36:14	1.16	217	2.50	-	-	*	
Total Hepta-Dioxins	1.77e+07			42:55	0.99	108	2.50	-	-	8	
Total Tetra-Furans	2.13e+07			23:10	1.34	40.6	2.50	-	-	*	
1st Fn. Tot Penta-Furans	1.71e+07			28:32	0.92	51.0	2.50	-	-	*	
Total Penta-Furans	4.85e+07			30:18	0.92	145	2.50	-	-	PeCDF 1	
Total Hexa-Furans	8.03e+07			35:21	1.34	277	2.50	-	-	*	
Total Hepta-Furans	2.00e+07			42:24	1.06	107	2.50	-	-	196 13	
										12	
										11	

Analyst: JC

Date: 3/6/09

FAL ID: ST022409M4      Filename: 24FEB09M      Sam:18      Acquired: 24-FEB-09 05:17:02      iCal: PCDDFAL3-2-5-09

Client ID: 1613 CS3 (080827J)      ConCal: ST022409M3      EndCal: ST022409M4

Results: 5340      GC Column: DB5      Amount: 1.000      NATO 1989 Tox: 100      WHO 1998 Tox: 126      WHO 2005 Tox: 115

Name	Resp	RA	RT	RRF	Conc	Qual	Fac	Noise-1	Noise-2	DL
------	------	----	----	-----	------	------	-----	---------	---------	----

2,3,7,8-TCDD	1.93e+06	0.81	y	27:31	1.04	9.78	2.50	-	-	*
1,2,3,7,8-PeCDD	8.24e+06	1.56	y	33:20	0.90	51.3	2.50	-	-	*
1,2,3,4,7,8-HxCDD	9.13e+06	1.27	y	38:41	1.43	48.6	2.50	-	-	*
1,2,3,6,7,8-HxCDD	6.88e+06	1.29	y	38:51	1.02	47.8	2.50	-	-	*
1,2,3,7,8,9-HxCDD	7.48e+06	1.29	y	39:18	1.05	52.2	2.50	-	-	*
1,2,3,4,6,7,8-HpCDD	5.28e+06	1.04	y	44:17	0.99	45.2	2.50	-	-	*
OCDD	7.18e+06	0.88	y	49:52	0.94	93.5	2.50	-	-	*

2,3,7,8-TCDF	3.73e+06	0.71	y	26:46	1.34	8.79	2.50	-	-	*
1,2,3,7,8-PeCDF	1.44e+07	1.62	y	31:36	1.03	49.7	2.50	-	-	*
2,3,4,7,8-PeCDF	1.11e+07	1.60	y	32:55	0.81	47.6	2.50	-	-	*
1,2,3,4,7,8-HxCDF	1.24e+07	1.26	y	37:18	1.35	53.1	2.50	-	-	*
1,2,3,6,7,8-HxCDF	1.17e+07	1.25	y	37:29	1.25	52.4	2.50	-	-	*
2,3,4,6,7,8-HxCDF	1.29e+07	1.25	y	38:26	1.53	51.0	2.50	-	-	*
1,2,3,7,8,9-HxCDF	8.41e+06	1.28	y	39:53	1.19	51.6	2.50	-	-	*
1,2,3,4,6,7,8-HpCDF	8.17e+06	1.02	y	42:22	1.10	52.7	2.50	-	-	*
1,2,3,4,7,8,9-HpCDF	6.50e+06	1.03	y	45:13	1.02	54.0	2.50	-	-	*
OCDF	8.25e+06	0.91	y	50:15	0.76	97.8	2.50	-	-	*

Rec

13C-2,3,7,8-TCDD	1.90e+07	0.78	y	27:30	0.94	103		103
13C-1,2,3,7,8-PeCDD	1.78e+07	1.77	y	33:18	0.75	120		120
13C-1,2,3,4,7,8-HxCDD	1.31e+07	1.29	y	38:40	1.03	92.7		92.7
13C-1,2,3,6,7,8-HxCDD	1.41e+07	1.28	y	38:50	1.08	94.9		94.9
13C-1,2,3,4,6,7,8-HpCDD	1.18e+07	1.06	y	44:15	0.89	96.1		96.1
13C-OCDD	1.63e+07	0.93	y	49:50	0.64	186		93.2

13C-2,3,7,8-TCDF	3.17e+07	0.85	y	26:45	0.86	103		103
13C-1,2,3,7,8-PeCDF	2.82e+07	1.65	y	31:34	0.76	104		104
13C-2,3,4,7,8-PeCDF	2.89e+07	1.66	y	32:53	0.78	104		104
13C-1,2,3,4,7,8-HxCDF	1.72e+07	0.53	y	37:16	1.34	93.0		93.0
13C-1,2,3,6,7,8-HxCDF	1.79e+07	0.54	y	37:29	1.40	92.5		92.5
13C-2,3,4,6,7,8-HxCDF	1.65e+07	0.53	y	38:24	1.29	93.1		93.1
13C-1,2,3,7,8,9-HxCDF	1.37e+07	0.54	y	39:51	1.01	99.0		99.0
13C-1,2,3,4,6,7,8-HpCDF	1.41e+07	0.44	y	42:21	1.16	88.4		88.4
13C-1,2,3,4,7,8,9-HpCDF	1.18e+07	0.45	y	45:11	0.96	89.1		89.1
13C-OCDF	2.21e+07	0.92	y	50:13	0.95	169		84.6

37Cl-2,3,7,8-TCDD	1.73e+06			27:31	0.88	9.97		99.7
-------------------	----------	--	--	-------	------	------	--	------

13C-1,2,3,4-TCDD	1.97e+07	0.79	y	26:54	-	53.1					
13C-1,2,3,4-TCDF	3.57e+07	0.84	y	25:39	-	59.5					
13C-1,2,3,7,8,9-HxCDD	1.38e+07	1.28	y	39:16	-	68.6					
Total Tetra-Dioxins	1.04e+07			23:31	1.04	52.7	2.50	-	-	*	20
Total Penta-Dioxins	2.74e+07			30:21	0.90	171	2.50	-	-	*	8
Total Hexa-Dioxins	3.49e+07			36:13	1.16	220	2.50	-	-	*	9
Total Hepta-Dioxins	1.28e+07			42:21	0.99	109	2.50	-	-	*	9
Total Tetra-Furans	1.71e+07			23:09	1.34	40.4	2.50	-	-	*	20
1st Fn. Tot Penta-Furans	1.43e+07			28:31	0.92	54.6	2.50	-	-	*	PeCDF 1
Total Penta-Furans	3.77e+07			30:17	0.92	144	2.50	-	-	*	198 9
Total Hexa-Furans	6.07e+07			35:20	1.34	278	2.50	-	-	*	7
Total Hepta-Furans	1.48e+07			42:22	1.06	108	2.50	-	-	*	12

Analyst: JK

Date: 3/2/09

FAL ID: ST030409A1      Filename: 04MAR09A      Sam:1      Acquired: 4-MAR-09 09:46:39      iCal: TCDFFAL1-1-13-09  
Client ID: 1613 CS3 (080827J)      ConCal: ST030409A1 EndCal: ST030409A2  
Results:      GC Column: DB225      Amount: 1.000

Name	Resp	RA	RT	RRF	Conc	Qual	Fac	Noise	DL	#Hom	Rec
2,3,7,8-TCDF	1.49e+07	0.78	y	19:34	1.03	8.94	2.50	-	-	1	104
13C-2,3,7,8-TCDF	1.62e+08	0.79	y	19:33	0.87	104					
13C-1,2,3,4-TCDF	1.79e+08	0.79	y	16:59	-	96.0					

Analyst: J Date: 3/4/09

FAL ID: ST030409A2      Filename: 04MAR09A      Sam:5      Acquired: 4-MAR-09 12:09:12      iCal: TCDFFAL1-1-13-09  
Client ID: 1613 CS3 (080827J)      ConCal: ST030409A1 EndCal: ST030409A2  
Results: 5340TCDF      GC Column: DB225      Amount: 1.000

Name	Resp	RA	RT	RRF	Conc	Qual	Fac	Noise	DL	#Hom
2,3,7,8-TCDF	9.51e+06	0.79	y	19:29	1.03	8.20	2.50	-	-	1

13C-2,3,7,8-TCDF	1.13e+08	0.78	y	19:27	0.87	93.1		Rec		93.1
13C-1,2,3,4-TCDF	1.39e+08	0.78	y	16:54	-	74.9				

Analyst: J Date: 3/4/09

1DFA - FORM I-HR CDD-1  
 CDD/CDF SAMPLE DATA SUMMARY  
 HIGH RESOLUTION

SAMPLE No.  
 DFBLK01

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564  
 LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: SDG NO.: SDG E3YC9  
 MATRIX: (SOIL/WATER/ASH/TISSUE/OIL) SOIL LAB SAMPLE ID: 1652-001-0001-MB  
 SAMPLE wt/vol: 1.000 (g/mL): g LAB FILE ID: 24FEB09M Sam: 13  
 WATER SAMPLE PREP: SPE (SEPF/SPE)  
 DATE RECEIVED: 18-FEB-09  
 CONCENTRATED EXTRACT VOLUME: 20 (uL)  
 DATE EXTRACTED: 18-FEB-09  
 INJECTION VOLUME: 2 (uL) % SOLIDS/LIPIDS: 100  
 DATE ANALYZED: 24-FEB-09  
 GC COLUMN: DB5 ID: 0.25 (mm)  
 DILUTION FACTOR: 1  
 CONCENTRATION UNITS: (pg/L or ng/kg) ng/kg

TARGET ANALYTE	SELECTED IONS	PEAK RT	ION RATIO #	CONCENTRATION	Q	EMPC/EDL
2,3,7,8-TCDD	320/322	NotFnd	*	*	U	0.605
2,3,7,8-TCDF	304/306	NotFnd	*	*	U	0.483
1,2,3,7,8-PeCDF	340/342	NotFnd	*	*	U	1.05
1,2,3,7,8-PeCDD	356/358	NotFnd	*	*	U	1.15
2,3,4,7,8-PeCDF	340/342	NotFnd	*	*	U	1.41
1,2,3,4,7,8-HxCDF	374/376	NotFnd	*	*	U	0.577
1,2,3,6,7,8-HxCDF	374/376	NotFnd	*	*	U	0.626
1,2,3,4,7,8-HxCDD	390/392	NotFnd	*	*	U	0.964
1,2,3,6,7,8-HxCDD	390/392	NotFnd	*	*	U	1.30
1,2,3,7,8,9-HxCDD	390/392	NotFnd	*	*	U	1.28
2,3,4,6,7,8-HxCDF	374/376	NotFnd	*	*	U	0.533
1,2,3,7,8,9-HxCDF	374/376	NotFnd	*	*	U	0.838
1,2,3,4,6,7,8-HpCDF	408/410	NotFnd	*	*	U	0.817
1,2,3,4,6,7,8-HpCDD	424/426	NotFnd	*	*	U	1.61
1,2,3,4,7,8,9-HpCDF	408/410	NotFnd	*	*	U	0.869
OCDD	458/460	NotFnd	*	*	U	3.93
OCDF	442/444	NotFnd	*	*	U	1.90

NOTE: Concentrations, Estimated Maximum Possible Concentrations (EMPCs), and Estimated Detection Levels (EDLs) for solid samples are calculated on a dry weight basis (except tissues, which are reported on a wet weight basis with % Lipids).

LABELED COMPOUNDS	SELECTED IONS	PEAK RT	ION RATIO #	ION RATIO LIMITS	% REC #	RECOVERY LIMITS
13C-2,3,7,8-TCDD	332/334	27:31	0.78	0.65-0.89	91.0	25-164
13C-1,2,3,7,8-PeCDD	368/370	33:19	1.77	1.32-1.78	101	25-181
13C-1,2,3,4,7,8-HxCDD	402/404	38:40	1.30	1.05-1.43	84.7	32-141
13C-1,2,3,6,7,8-HxCDD	402/404	38:51	1.29	1.05-1.43	83.5	28-130
13C-1,2,3,4,6,7,8-HpCDD	436/438	44:16	1.08	0.88-1.20	88.0	23-140
13C-OCDD	470/472	49:51	0.93	0.76-1.02	72.7	17-157
13C-2,3,7,8-TCDF	316/318	26:46	0.85	0.65-0.89	91.9	24-169
13C-1,2,3,7,8-PeCDF	352/354	31:35	1.66	1.32-1.78	91.4	24-185
13C-2,3,4,7,8-PeCDF	352/354	32:54	1.66	1.32-1.78	90.0	21-178
13C-1,2,3,4,7,8-HxCDF	384/386	37:17	0.53	0.43-0.59	80.0	26-152
13C-1,2,3,6,7,8-HxCDF	384/386	37:28	0.53	0.43-0.59	78.0	26-123
13C-1,2,3,7,8,9-HxCDF	384/386	39:52	0.54	0.43-0.59	88.4	29-147
13C-2,3,4,6,7,8-HxCDF	384/386	38:25	0.54	0.43-0.59	81.0	28-136
13C-1,2,3,4,6,7,8-HpCDF	418/420	42:22	0.44	0.37-0.51	74.8	28-143
13C-1,2,3,4,7,8,9-HpCDF	418/420	45:12	0.45	0.37-0.51	97.1	26-138
13C-OCDF	454/456	50:13	0.92	0.76-1.02	72.9	17-157
37Cl-2,3,7,8-TCDD	328/NA	27:32	NA	NA	94.1	35-197

# Column to be used to flag values outside (QC) limits.

ANALYST: E

DATE: 3/2/09

1DFB - FORM I-HR CDD-2  
 CDD/CDF TOXICITY EQUIVALENCE SUMMARY  
 HIGH RESOLUTION

SAMPLE No.  
 DFBLK01

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: SDG NO.: SDG E3YC9

MATRIX: (SOIL/WATER/ASH/TISSUE/OIL) SOIL LAB SAMPLE ID: 1652-001-0001-MB

SAMPLE wt/vol: 1.000 (g/mL): g LAB FILE ID: 24FEB09M Sam: 13

WATER SAMPLE PREP: SPE (SEPF/SPE) DATE RECEIVED: 18-FEB-09

CONCENTRATED EXTRACT VOLUME: 20 (uL) DATE EXTRACTED: 18-FEB-09

INJECTION VOLUME: 2 (uL) % SOLIDS/LIPIDS: 100 DATE ANALYZED: 24-FEB-09

GC COLUMN: DB5 ID: 0.25 (mm) DILUTION FACTOR: 1

CONCENTRATION UNITS: (pg/L or ng/kg) ng/kg


TARGET ANALYTE	CONCENTRATION	TEF*	TEF-ADJUSTED CONCENTRATION
2,3,7,8-TCDD	*	x 1.0 =	*
2,3,7,8-TCDF	*	x 0.1 =	*
1,2,3,7,8-PeCDD	*	x 0.05 =	*
1,2,3,7,8-PeCDD	*	x 0.5 =	*
2,3,4,7,8-PeCDF	*	x 0.5 =	*
1,2,3,4,7,8-HxCDF	*	x 0.1 =	*
1,2,3,6,7,8-HxCDF	*	x 0.1 =	*
1,2,3,4,7,8-HxCDD	*	x 0.1 =	*
1,2,3,6,7,8-HxCDD	*	x 0.1 =	*
1,2,3,7,8,9-HxCDD	*	x 0.1 =	*
2,3,4,6,7,8-HxCDF	*	x 0.1 =	*
1,2,3,7,8,9-HxCDF	*	x 0.1 =	*
1,2,3,4,6,7,8-HpCDF	*	x 0.01 =	*
1,2,3,4,6,7,8-HpCDD	*	x 0.01 =	*
1,2,3,4,7,8,9-HpCDF	*	x 0.01 =	*
OCDD	*	x 0.001 =	*
OCDF	*	x 0.001 =	*
Total =		0.00	

\* TEF - Toxicity Equivalent Factors from EPA/625/3-89/016 March 1989 - Interim Procedures for Estimating Risks Associated with Exposures to Mixtures of Chlorinated Dibenzo-p-Dioxins and Chlorinated Dibenzofurans (CDDs and CDFs) and 1989 Update


ANALYST: E

DATE: 3/2/09


File:24FEB09M #1-390 Acq:24-FEB-2009 00:40:37 GC EI+ Voltage SIR Autospec-Ultima  
319.8965 S:13 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-390 Acq:24-FEB-2009 00:40:37 GC EI+ Voltage SIR Autospec-Ultima  
321.8936 S:13 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-390 Acq:24-FEB-2009 00:40:37 GC EI+ Voltage SIR Autospec-Ultima  
331.9368 S:13 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-390 Acq:24-FEB-2009 00:40:37 GC EI+ Voltage SIR Autospec-Ultima  
333.9339 S:13 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-425 Acq:24-FEB-2009 00:40:37 GC EI+ Voltage SIR Autospec-Ultima  
355.8546 S:13 R:2 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100.0,0.00%,FF) Exp:PCDD  
Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TG



File:24FEB09M #1-425 Acq:24-FEB-2009 00:40:37 GC EI+ Voltage SIR Autospec-Ultima  
357.8517 S:13 R:2 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100.0,0.00%,FF) Exp:PCDD  
Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TG







File:24FEB09M #1-425 Acq:24-FEB-2009 00:40:37 GC EI+ Voltage SIR Autospec-Ultima  
367.8949 S:13 R:2 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100.0,0.00%,FF) Exp:PCDD  
Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TG







File:24FEB09M #1-425 Acq:24-FEB-2009 00:40:37 GC EI+ Voltage SIR Autospec-Ultima  
LOCK MASS CHECK S:13 R:2 Exp:PCDD  
Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TG

File:24FEB09M #1-464 Acq:24-FEB-2009 00:40:37 GC EI+ Voltage SIR Autospec-Ultima  
389.8156 S:13 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100,0,0.00%,FF) Exp:PCDD





Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-541 Acq:24-FEB-2009 00:40:37 GC EI+ Voltage SIR Autospec-Ultima  
423.7767 S:13 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-347 Acq:24-FEB-2009 00:40:37 GC EI+ Voltage SIR Autospec-Ultima  
457.7377 S:13 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-390 Acq:24-FEB-2009 00:40:37 GC EI+ Voltage SIR Autospec-Ultima  
 303.9016 S:13 BSUB(10000,15,-3.0) PKD(5,5,3,0 10%,100.0,0.00%,F,F) Exp:PCDD  
 Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-390 Acq:24-FEB-2009 00:40:37 GC EI+ Voltage SIR Autospec-Ultima  
 305.8987 S:13 BSUB(10000,15,-3.0) PKD(5,5,3,0 10%,100.0,0.00%,F,F) Exp:PCDD  
 Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TC

File:24FEB09M #1-390 Acq:24-FEB-2009 00:40:37 GC EI+ Voltage SIR Autospec-Ultima  
 315.9419 S:13 BSUB(10000,15,-3.0) PKD(5,5,3,0 10%,100.0,0.00%,F,F) Exp:PCDD  
 Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TC

File:24FEB09M #1-390 Acq:24-FEB-2009 00:40:37 GC EI+ Voltage SIR Autospec-Ultima  
 317.9389 S:13 BSUB(10000,15,-3.0) PKD(5,5,3,0 10%,100.0,0.00%,F,F) Exp:PCDD  
 Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-390 Acq:24-FEB-2009 00:40:37 GC EI+ Voltage SIR Autospec-Ultima  
 375.8364 S:13 BSUB(10000,15,-3.0) PKD(5,5,3,0 10%,100.0,0.00%,F,F) Exp:PCDD  
 Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-390 Acq:24-FEB-2009 00:40:37 GC EI+ Voltage SIR Autospec-Ultima  
 375.8364 S:13 BSUB(10000,15,-3.0) PKD(5,5,3,0 10%,100.0,0.00%,F,F) Exp:PCDD  
 Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TC

File:24FEB09M #1-390 Acq:24-FEB-2009 00:40:37 GC EI+ Voltage SIR Autospec-Ultima  
 375.8364 S:13 BSUB(10000,15,-3.0) PKD(5,5,3,0 10%,100.0,0.00%,F,F) Exp:PCDD  
 Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TC


File:24FEB09M #1-390 Acq:24-FEB-2009 00:40:37 GC EI + Voltage SIR Autospec-Ultima  
339.8597 S:13 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-390 Acq:24-FEB-2009 00:40:37 GC EI + Voltage SIR Autospec-Ultima  
341.8568 S:13 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-390 Acq:24-FEB-2009 00:40:37 GC EI + Voltage SIR Autospec-Ultima  
409.7974 S:13 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-390 Acq:24-FEB-2009 00:40:37 GC EI + Voltage SIR Autospec-Ultima  
330.9792 S:13 Exp:PCDD  
Sample Text:Frontier Analytical Laboratory Analyst ID:TC


File:24FEB09M #1-425 Acq:24-FEB-2009 00:40:37 GC EI+ Voltage SIR Autospec-Ultima  
 339.8597 S:13 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
 Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-425 Acq:24-FEB-2009 00:40:37 GC EI+ Voltage SIR Autospec-Ultima  
 341.8568 S:13 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
 Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-425 Acq:24-FEB-2009 00:40:37 GC EI+ Voltage SIR Autospec-Ultima  
 353.8970 S:13 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
 Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-425 Acq:24-FEB-2009 00:40:37 GC EI+ Voltage SIR Autospec-Ultima  
 409.7974 S:13 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
 Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-464 Acq:24-FEB-2009 00:40:37 GC EI+ Voltage SIR Autospec-Ultima  
373.8207 S:13 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-464 Acq:24-FEB-2009 00:40:37 GC EI+ Voltage SIR Autospec-Ultima  
383.8639 S:13 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-464 Acq:24-FEB-2009 00:40:37 GC EI+ Voltage SIR Autospec-Ultima  
383.8639 S:13 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-464 Acq:24-FEB-2009 00:40:37 GC EI+ Voltage SIR Autospec-Ultima  
445.7555 S:13 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-541 Acq:24-FEB-2009 00:40:37 GC EI+ Voltage SIR Autospec-Ultima  
407.7818 S:13 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD

Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TC



Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TC



Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TC

42:22  
A1.14E7  
A5.05E6


42:24  
A3.14E3  
A1.42E3  
A911.47  
A1.04E3  
A1.25E3  
A4.40E3  
A1.19E3  
A1.41E3  
A880.97  
A537.46  
A695.80  
A5.47E6

File:24FEB09M #1-541 Acq:24-FEB-2009 00:40:37 GC EI+ Voltage SIR Autospec-Ultima  
419.8220 S:13 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD

Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TC

42:22  
A1.22E7  
A45.12


A1.22E7  
A45.12




Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TC

42:22  
A2.22E3  
A862.11  
A1.15E3  
A983.37  
A798.28  
A1.29E3  
A1.38E3  
A3.18E3  
A1.39E3  
1.4E3  
7.0E2  
A687.29  
0.0E0


File:24FEB09M #1-347 Acq:24-FEB-2009 00:40:37 GC EI+ Voltage SIR Autospec-Ultima  
441.7428 S:13 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-347 Acq:24-FEB-2009 00:40:37 GC EI+ Voltage SIR Autospec-Ultima  
443.7398 S:13 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-347 Acq:24-FEB-2009 00:40:37 GC EI+ Voltage SIR Autospec-Ultima  
453.7831 S:13 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-347 Acq:24-FEB-2009 00:40:37 GC EI+ Voltage SIR Autospec-Ultima  
513.6775 S:13 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:1652-001-0001-MB File Text:Frontier Analytical Laboratory Analyst ID:TC



FAL ID: 1652-001-0001-MB      Filename: 24FEB09M      Sam:13      Acquired: 24-FEB-09 00:40:37      iCal: PCDDFAL3-2-5-09  
 Client ID: DFBLK01      ConCal: ST022409M3      EndCal: ST022409M4  
 Results: 5340      GC Column: db5      Amount: 1.000      NATO 1989 Tox: 0.00  
 WHO 1998 Tox: 0.00      WHO 2005 Tox: 0.00  
 Name      Resp      RA      RT      RRF      Conc      Qual      Fac Noise-1      Noise-2      DL

2,3,7,8-TCDD	*	* n NotFnd	1.04	*	U	2.50	207	454	0.605	
1,2,3,7,8-PeCDD	*	* n NotFnd	0.90	*	U	2.50	647	213	1.15	
1,2,3,4,7,8-HxCDD	*	* n NotFnd	1.43	*	U	2.50	547	357	0.964	
1,2,3,6,7,8-HxCDD	*	* n NotFnd	1.02	*	U	2.50	547	357	1.30	
1,2,3,7,8,9-HxCDD	*	* n NotFnd	1.05	*	U	2.50	547	357	1.28	
1,2,3,4,6,7,8-HpCDD	*	* n NotFnd	0.99	*	U	2.50	510	383	1.61	
OCDD	*	* n NotFnd	0.94	*	U	2.50	443	511	3.93	
2,3,7,8-TCDF	*	* n NotFnd	1.34	*	U	2.50	403	835	0.483	
1,2,3,7,8-PeCDF	*	* n NotFnd	1.03	*	U	2.50	412	1170	1.05	
2,3,4,7,8-PeCDF	*	* n NotFnd	0.81	*	U	2.50	412	1170	1.41	
1,2,3,4,7,8-HxCDF	*	* n NotFnd	1.35	*	U	2.50	335	294	0.577	
1,2,3,6,7,8-HxCDF	*	* n NotFnd	1.25	*	U	2.50	335	294	0.626	
2,3,4,6,7,8-HxCDF	*	* n NotFnd	1.53	*	U	2.50	335	294	0.533	
1,2,3,7,8,9-HxCDF	*	* n NotFnd	1.19	*	U	2.50	335	294	0.838	
1,2,3,4,6,7,8-HpCDF	*	* n NotFnd	1.10	*	U	2.50	244	333	0.817	
1,2,3,4,7,8,9-HpCDF	*	* n NotFnd	1.02	*	U	2.50	244	333	0.869	
OCDF	*	* n NotFnd	0.76	*	U	2.50	225	336	1.90	
								Rec		
13C-2,3,7,8-TCDD	2.23e+07	0.78 y	27:31	0.94	1820				91.0	
13C-1,2,3,7,8-PeCDD	1.99e+07	1.77 y	33:19	0.75	2020				101	
13C-1,2,3,4,7,8-HxCDD	1.66e+07	1.30 y	38:40	1.03	1690				84.7	
13C-1,2,3,6,7,8-HxCDD	1.71e+07	1.29 y	38:51	1.08	1670				83.5	
13C-1,2,3,4,6,7,8-HpCDD	1.50e+07	1.08 y	44:16	0.89	1760				88.0	
13C-OCDD	1.76e+07	0.93 y	49:51	0.64	2910				72.7	
13C-2,3,7,8-TCDF	3.75e+07	0.85 y	26:46	0.86	1840				91.9	
13C-1,2,3,7,8-PeCDF	3.28e+07	1.66 y	31:35	0.76	1830				91.4	
13C-2,3,4,7,8-PeCDF	3.30e+07	1.66 y	32:54	0.78	1800				90.0	
13C-1,2,3,4,7,8-HxCDF	2.04e+07	0.53 y	37:17	1.34	1600				80.0	
13C-1,2,3,6,7,8-HxCDF	2.08e+07	0.53 y	37:28	1.40	1560				78.0	
13C-2,3,4,6,7,8-HxCDF	1.99e+07	0.54 y	38:25	1.29	1620				81.0	
13C-1,2,3,7,8,9-HxCDF	1.69e+07	0.54 y	39:52	1.01	1770				88.4	
13C-1,2,3,4,6,7,8-HpCDF	1.65e+07	0.44 y	42:22	1.16	1500				74.8	
13C-1,2,3,4,7,8,9-HpCDF	1.77e+07	0.45 y	45:12	0.96	1940				97.1	
13C-OCDF	2.63e+07	0.92 y	50:13	0.95	2920				72.9	
37Cl-2,3,7,8-TCDD	8.64e+06		27:32	0.88	753				94.1	
13C-1,2,3,4-TCDD	2.61e+07	0.79 y	26:55	-	70.3					
13C-1,2,3,4-TCDF	4.73e+07	0.86 y	25:40	-	78.9					
13C-1,2,3,7,8,9-HxCDD	1.90e+07	1.29 y	39:17	-	94.7					
Total Tetra-Dioxins	*	NotFnd	1.04	*	U	2.50	207	454	0.605	0
Total Penta-Dioxins	*	NotFnd	0.90	*	U	2.50	647	213	1.15	0
Total Hexa-Dioxins	*	NotFnd	1.16	*	U	2.50	547	357	1.30	0
Total Hepta-Dioxins	*	NotFnd	0.99	*	U	2.50	510	383	1.61	0
Total Tetra-Furans	*	NotFnd	1.34	*	U	2.50	403	835	0.483	0
1st Fn. Tot Penta-Furans	*	NotFnd	0.92	*	U	2.50	412	1170	1.41	PeCDF 0
Total Penta-Furans	*	NotFnd	0.92	*	U	2.50	412	1170	1.41	0.00 0
Total Hexa-Furans	*	NotFnd	1.34	*	U	2.50	335	294	0.838	0
Total Hepta-Furans	*	NotFnd	1.06	*	U	2.50	244	333	0.869	0

Analyst: SC

Date: 3/2/09

2DF - FORM II-HR CDD  
CDD/CDF TOTAL HOMOLOGUE CONCENTRATION SUMMARY  
HIGH RESOLUTION

SAMPLE No.  
DFBLK01

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564

LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: SDG NO.: SDG E3YC9

MATRIX: (SOIL/WATER/ASH/TISSUE/OIL) SOIL LAB SAMPLE ID: 1652-001-0001-MB

SAMPLE wt/vol: 1.000 (g/mL): g LAB FILE ID: 24FEB09M Sam: 13

WATER SAMPLE PREP: SPE (SEPF/SPE) DATE RECEIVED: 18-FEB-09

CONCENTRATED EXTRACT VOLUME: 20 (uL) DATE EXTRACTED: 18-FEB-09

INJECTION VOLUME: 2 (uL) % SOLIDS/LIPIDS: 100 DATE ANALYZED: 24-FEB-09

GC COLUMN: DB5 ID: 0.25 (mm) DILUTION FACTOR: 1

CONCENTRATION UNITS: (pg/L or ng/kg) ng/kg

HOMOLOGUE	PEAKS	CONCENTRATION	Q	EMPC/EDL
<b>DIOXINS</b>				
Total TCDD	0	*	U	0.605
Total PeCDD	0	*	U	1.15
Total HxCDD	0	*	U	1.30
Total HpCDD	0	*	U	1.61
<b>FURANS</b>				
Total TCDF	0	*	U	0.483
Total PeCDF	0	*	U	1.41
Total HxCDF	0	*	U	0.838
Total HpCDF	0	*	U	0.869

NOTE: Concentrations, Estimated Maximum Possible Concentrations (EMPCs), and Estimated Detection Limits (EDLs) for solid samples are calculated on a dry weight basis (except tissues, which are reported on a wet weight basis with % Lipids). The total homologue concentrations do not affect the TEF (Toxicity Equivalent Factor) calculations.

ANALYST: SC

DATE: 3/2/09

1DFA - FORM I-HR CDD-1  
 CDD/CDF SAMPLE DATA SUMMARY  
 HIGH RESOLUTION

SAMPLE No.  
 DLCS01

LAB NAME: FRONTIER ANALYTICAL LAB CONTRACT: EP08W001564  
 LAB CODE: FALE CASE NO.: EPA05-2008-01 TO NO.: SDG NO.: SDG E3YC9  
 MATRIX: (SOIL/WATER/ASH/TISSUE/OIL) SOIL LAB SAMPLE ID: 1652-001-0001-OPR  
 SAMPLE wt/vol: 1.000 (g/mL): g LAB FILE ID: 24FEB09M Sam: 12  
 WATER SAMPLE PREP: SPE (SEPF/SPE)  
 DATE RECEIVED: 11-FEB-09  
 CONCENTRATED EXTRACT VOLUME: 20 (uL)  
 DATE EXTRACTED: 11-FEB-09  
 INJECTION VOLUME: 2 (uL) % SOLIDS/LIPIDS: 100  
 DATE ANALYZED: 23-FEB-09  
 GC COLUMN: DB5 ID: 0.25 (mm)  
 DILUTION FACTOR: 1  
 CONCENTRATION UNITS: (pg/L or ng/kg) ng/kg

TARGET ANALYTE	SELECTED IONS	PEAK RT	ION RATIO #	CONCENTRATION	Q	EMPC/EDL
2,3,7,8-TCDD	320/322	27:32	0.79	9.67	*	
2,3,7,8-TCDF	304/306	26:47	0.70	8.38	*	
1,2,3,7,8-PeCDD	340/342	31:36	1.60	46.1	*	
1,2,3,7,8-PeCDD	356/358	33:20	1.59	50.1	*	
2,3,4,7,8-PeCDF	340/342	32:55	1.60	46.4	*	
1,2,3,4,7,8-HxCDF	374/376	37:18	1.27	51.1	*	
1,2,3,6,7,8-HxCDF	374/376	37:31	1.26	50.5	*	
1,2,3,4,7,8-HxCDD	390/392	38:41	1.28	47.0	*	
1,2,3,6,7,8-HxCDD	390/392	38:51	1.27	46.0	*	
1,2,3,7,8,9-HxCDD	390/392	39:18	1.27	46.2	*	
2,3,4,6,7,8-HxCDF	374/376	38:26	1.27	50.1	*	
1,2,3,7,8,9-HxCDF	374/376	39:53	1.28	50.4	*	
1,2,3,4,6,7,8-HpCDF	408/410	42:23	1.02	50.8	*	
1,2,3,4,6,7,8-HpCDD	424/426	44:17	1.06	43.4	*	
1,2,3,4,7,8,9-HpCDF	408/410	45:13	1.03	51.8	*	
OCDD	458/460	49:53	0.90	89.2	*	
OCDF	442/444	50:15	0.91	94.3	*	

NOTE: Concentrations, Estimated Maximum Possible Concentrations (EMPCs), and Estimated Detection Levels (EDLs) for solid samples are calculated on a dry weight basis (except tissues, which are reported on a wet weight basis with % Lipids).





LABELED COMPOUNDS	SELECTED IONS	PEAK RT	ION RATIO #	ION RATIO LIMITS	% REC #	RECOVERY LIMITS
13C-2,3,7,8-TCDD	332/334	27:31	0.78	0.65-0.89	100	25-164
13C-1,2,3,7,8-PeCDD	368/370	33:19	1.78	1.32-1.78	114	25-181
13C-1,2,3,4,7,8-HxCDD	402/404	38:40	1.29	1.05-1.43	93.2	32-141
13C-1,2,3,6,7,8-HxCDD	402/404	38:50	1.28	1.05-1.43	93.0	28-130
13C-1,2,3,4,6,7,8-HpCDD	436/438	44:17	1.07	0.88-1.20	96.9	23-140
13C-OCDD	470/472	49:51	0.93	0.76-1.02	81.2	17-157
13C-2,3,7,8-TCDF	316/318	26:45	0.86	0.65-0.89	101	24-169
13C-1,2,3,7,8-PeCDF	352/354	31:35	1.66	1.32-1.78	103	24-185
13C-2,3,4,7,8-PeCDF	352/354	32:54	1.65	1.32-1.78	101	21-178
13C-1,2,3,4,7,8-HxCDF	384/386	37:17	0.53	0.43-0.59	90.2	26-152
13C-1,2,3,6,7,8-HxCDF	384/386	37:29	0.53	0.43-0.59	88.9	26-123
13C-1,2,3,7,8,9-HxCDF	384/386	39:52	0.53	0.43-0.59	97.2	29-147
13C-2,3,4,6,7,8-HxCDF	384/386	38:26	0.53	0.43-0.59	89.5	28-136
13C-1,2,3,4,6,7,8-HpCDF	418/420	42:22	0.44	0.37-0.51	83.5	28-143
13C-1,2,3,4,7,8,9-HpCDF	418/420	45:12	0.44	0.37-0.51	108	26-138
13C-OCDF	454/456	50:13	0.92	0.76-1.02	82.3	17-157
37Cl-2,3,7,8-TCDD	328/NA	27:32	NA	NA	101	35-197

# Column to be used to flag values outside (QC) limits.




ANALYST: 

DATE:  3/1/09




File:24FEB09M #1-390 Acq:23-FEB-2009 23:45:23 GC EI+ Voltage SIR Autospec-Ultima  
319.8965 S:12 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:1652-001-0001-OPR File Text:Frontier Analytical Laboratory Analyst ID:TC







File:24FEB09M #1-425 Acq:23-FEB-2009 23:45:23 GC EI+ Voltage SIR Autospec-Ultima  
355.8546 S:12 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:1652-001-0001-OPR File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-464 Acq:23-FEB-2009 23:45:23 GC EI+ Voltage SIR Autospec-Ultima  
389.8156 S:12 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:1652-001-0001-OPR File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-541 Acq:23-FEB-2009 23:45:23 GC EI+ Voltage SIR Autospec-Ultima  
423.7767 S:12 F:4 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:1652-001-0001-OPR File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-347 Acq:23-FEB-2009 23:45:23 GC EI+ Voltage SIR Autospec-Ultima  
457.7377 S:12 F:5 BSUB10000,15,-3.0 PKD(5,5,3,0.10%,100,0,0.00%,FF) Exp:PCDD  
Sample Text:1652-001-0001-OPR File Text:Frontier Analytical Laboratory Analyst ID:TC

100 %  
50  
0

49:00 50:00 51:00 52:00  
Time  
A4.30E6  
49:53

File:24FEB09M #1-347 Acq:23-FEB-2009 23:45:23 GC EI+ Voltage SIR Autospec-Ultima  
459.7348 S:12 F:5 BSUB10000,15,-3.0 PKD(5,5,3,0.10%,100,0,0.00%,FF) Exp:PCDD  
Sample Text:1652-001-0001-OPR File Text:Frontier Analytical Laboratory Analyst ID:TC

100 %  
50  
0

49:00 50:00 51:00 52:00  
Time  
A4.80E6  
49:53

File:24FEB09M #1-347 Acq:23-FEB-2009 23:45:23 GC EI+ Voltage SIR Autospec-Ultima  
469.7780 S:12 F:5 BSUB10000,15,-3.0 PKD(5,5,3,0.10%,100,0,0.00%,FF) Exp:PCDD  
Sample Text:1652-001-0001-OPR File Text:Frontier Analytical Laboratory Analyst ID:TC

100 %  
50  
0

49:00 50:00 51:00 52:00  
Time  
A1.05E7  
49:51

File:24FEB09M #1-347 Acq:23-FEB-2009 23:45:23 GC EI+ Voltage SIR Autospec-Ultima  
471.7750 S:12 F:5 BSUB10000,15,-3.0 PKD(5,5,3,0.10%,100,0,0.00%,FF) Exp:PCDD  
Sample Text:1652-001-0001-OPR File Text:Frontier Analytical Laboratory Analyst ID:TC

100 %  
50  
0

49:00 50:00 51:00 52:00  
Time  
1.5E6  
49:51

File:24FEB09M #1-347 Acq:23-FEB-2009 23:45:23 GC EI+ Voltage SIR Autospec-Ultima  
474.9728 S:12 F:5 Exp:PCDD  
Sample Text:1652-001-0001-OPR File Text:Frontier Analytical Laboratory Analyst ID:TC





100 %  
50  
0

49:00 50:00 51:00 52:00  
Time  
A1.12E7  
49:51

1.6E6  
8.2E5  
0.0E0




Time

File:24FEB09M #1-390 Acq:23-FEB-2009 23:45:23 GC EI+ Voltage SIR Autospec-Ultima  
303.9016 S:12 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:1652-001-0001-OPR File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-390 Acq:23-FEB-2009 23:45:23 GC EI+ Voltage SIR Autospec-Ultima  
339.8597 S:12 BSUB(10000,15,-3.0) PKD(5.5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:1652-001-0001-OPR File Text:Frontier Analytical Laboratory Analyst ID:TC





000264 of 000279




Time

0.0E0  
1.4E7  
2.7E7


File:24FEB09M #1-425 Acq:23-FEB-2009 23:45:23 GC EI+ Voltage SIR Autospec-Ultima  
339.8597 S:12 F:2 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100,0,0,0.00%,F,F) Exp:PCDD  
Sample Text:1652-001-0001-OPR File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-464 Acq:23-FEB-2009 23:45:23 GC EI+ Voltage SIR Autospec-Ultima  
373.8207 S:12 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:1652-001-0001-OPR File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-464 Acq:23-FEB-2009 23:45:23 GC EI+ Voltage SIR Autospec-Ultima  
383.8639 S:12 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:1652-001-0001-OPR File Text:Frontier Analytical Laboratory Analyst ID:TC



File:24FEB09M #1-464 Acq:23-FEB-2009 23:45:23 GC EI+ Voltage SIR Autospec-Ultima  
385.8610 S:12 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:1652-001-0001-OPR File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-464 Acq:23-FEB-2009 23:45:23 GC EI+ Voltage SIR Autospec-Ultima  
445.7555 S:12 F:3 BSUB(10000,15,-3.0) PKD(5,5,3,0.10%,100,0,0.00%,F,F) Exp:PCDD  
Sample Text:1652-001-0001-OPR File Text:Frontier Analytical Laboratory Analyst ID:TC




File:24FEB09M #1-541 Acq:23-FEB-2009 23:45:23 GC EI+ Voltage SIR Autospec-Ultima  
407.7818 S:12 F:4 BSUB10000.15,-3.0 PKD(5,5,3,0.10%,100.0,0.00%,F,F) Exp:PCDD  
Sample Text:1652-001-0001-OPR File Text:Frontier Analytical Laboratory Analyst ID:TC


42:23  
A5.73E6  
45:13  
A5.82E6  
1.1E6  
5.5E5  
0.0E0 Time




42:23  
A5.61E6  
45:12  
A5.64E6  
1.1E6  
5.4E5  
0.0E0 Time




45:12  
A6.63E6  
1.2E6  
6.0E5  
0.0E0 Time



42:22  
A6.24E6  
45:12  
A1.41E7  
1.2E6  
6.0E5  
0.0E0 Time






42:22  
A1.41E7  
45:12  
A1.50E7  
1.2E6  
6.0E5  
0.0E0 Time




42:00  
42.07  
A1.04E3  
42:43  
A1.13E3  
43:13  
A945.23  
43:40  
A1.73E3  
44:10  
A1.27E3  
45:02  
A1.13E3  
45:19  
A578.14  
45:57  
A1.37E3  
46:24  
A1.06E3  
47:05  
A1.04E3  
47:40  
A1.64E3  
5.0E2  
0.0E0 Time

File:24FEB09M #1-347 Acq:23-FEB-2009 23:45:23 GC EI+ Voltage SIR Autospec-Ultima  
441,7428 S:12 F:5 BSUB(10000,15,-3.0) PKD(5,5,3,0,10%,100,0,0,00%,F,F) Exp:PCDD  
Sample Text:1652-001-0001-OPR File Text:Frontier Analytical Laboratory Analyst ID:TC



FAL ID: 1652-001-0001-OPR      Filename: 24FEB09M      Sam:12      Acquired: 23-FEB-09 23:45:23      ICAL: PCDDFAL3-2-5-09  
 Client ID: DLCS01      ConCal: ST022409M3      EndCal: ST022409M4  
 Results: 5340      GC Column: db5      Amount: 1.000      NATO 1989 Tox: 96.8  
 WHO 1998 Tox: 122      WHO 2005 Tox: 112  
 Name      Resp      RA      RT      RRF      Conc      Qual      Fac Noise-1      Noise-2      DL

2,3,7,8-TCDD	2.70e+06	0.79	y	27:32	1.04	9.67	2.50	-	-	*		
1,2,3,7,8-PeCDD	1.11e+07	1.59	y	33:20	0.90	50.1	2.50	-	-	*		
1,2,3,4,7,8-HxCDD	1.35e+07	1.28	y	38:41	1.43	47.0	2.50	-	-	*		
1,2,3,6,7,8-HxCDD	9.90e+06	1.27	y	38:51	1.02	46.0	2.50	-	-	*		
1,2,3,7,8,9-HxCDD	1.00e+07	1.27	y	39:18	1.05	46.2	2.50	-	-	*		
1,2,3,4,6,7,8-HpCDD	7.80e+06	1.06	y	44:17	0.99	43.4	2.50	-	-	*		
OCDD	9.10e+06	0.90	y	49:53	0.94	89.2	2.50	-	-	*		
2,3,7,8-TCDF	5.06e+06	0.70	y	26:47	1.34	8.38	2.50	-	-	*		
1,2,3,7,8-PeCDF	1.92e+07	1.60	y	31:36	1.03	46.1	2.50	-	-	*		
2,3,4,7,8-PeCDF	1.52e+07	1.60	y	32:55	0.81	46.4	2.50	-	-	*		
1,2,3,4,7,8-HxCDF	1.77e+07	1.27	y	37:18	1.35	51.1	2.50	-	-	*		
1,2,3,6,7,8-HxCDF	1.66e+07	1.26	y	37:31	1.25	50.5	2.50	-	-	*		
2,3,4,6,7,8-HxCDF	1.86e+07	1.27	y	38:26	1.53	50.1	2.50	-	-	*		
1,2,3,7,8,9-HxCDF	1.23e+07	1.28	y	39:53	1.19	50.4	2.50	-	-	*		
1,2,3,4,6,7,8-HpCDF	1.13e+07	1.02	y	42:23	1.10	50.8	2.50	-	-	*		
1,2,3,4,7,8,9-HpCDF	1.15e+07	1.03	y	45:13	1.02	51.8	2.50	-	-	*		
OCDF	1.18e+07	0.91	y	50:15	0.76	94.3	2.50	-	-	*		
									Rec			
13C-2,3,7,8-TCDD	2.69e+07	0.78	y	27:31	0.94	100				100		
13C-1,2,3,7,8-PeCDD	2.46e+07	1.78	y	33:19	0.75	114				114		
13C-1,2,3,4,7,8-HxCDD	2.02e+07	1.29	y	38:40	1.03	93.2				93.2		
13C-1,2,3,6,7,8-HxCDD	2.11e+07	1.28	y	38:50	1.08	93.0				93.0		
13C-1,2,3,4,6,7,8-HpCDD	1.82e+07	1.07	y	44:17	0.89	96.9				96.9		
13C-OCDD	2.17e+07	0.93	y	49:51	0.64	162				81.2		
13C-2,3,7,8-TCDF	4.51e+07	0.86	y	26:45	0.86	101				101		
13C-1,2,3,7,8-PeCDF	4.04e+07	1.66	y	31:35	0.76	103				103		
13C-2,3,4,7,8-PeCDF	4.05e+07	1.65	y	32:54	0.78	101				101		
13C-1,2,3,4,7,8-HxCDF	2.55e+07	0.53	y	37:17	1.34	90.2				90.2		
13C-1,2,3,6,7,8-HxCDF	2.62e+07	0.53	y	37:29	1.40	88.9				88.9		
13C-2,3,4,6,7,8-HxCDF	2.42e+07	0.53	y	38:26	1.29	89.5				89.5		
13C-1,2,3,7,8,9-HxCDF	2.06e+07	0.53	y	39:52	1.01	97.2				97.2		
13C-1,2,3,4,6,7,8-HpCDF	2.03e+07	0.44	y	42:22	1.16	83.5				83.5		
13C-1,2,3,4,7,8,9-HpCDF	2.17e+07	0.44	y	45:12	0.96	108				108		
13C-OCDF	3.28e+07	0.92	y	50:13	0.95	165				82.3		
37CL-2,3,7,8-TCDD	1.02e+07			27:32	0.88	40.4				101		
13C-1,2,3,4-TCDD	2.86e+07	0.78	y	26:56	-	76.9						
13C-1,2,3,4-TCDF	5.15e+07	0.86	y	25:41	-	86.0						
13C-1,2,3,7,8,9-HxCDD	2.10e+07	1.29	y	39:17	-	105						
Total Tetra-Dioxins	2.81e+06			23:08	1.04	10.0	2.50	-	-	*	21	
Total Penta-Dioxins	1.11e+07			33:20	0.90	50.2	2.50	-	-	*	7	
Total Hexa-Dioxins	3.35e+07			38:41	1.16	139	2.50	-	-	*	7	
Total Hepta-Dioxins	7.89e+06			42:55	0.99	43.9	2.50	-	-	*	8	
1st Fn. Tot Penta-Furans	5.25e+06			24:19	1.34	8.70	2.50	-	-	*	10	
Total Penta-Furans	7.59e+04			22:41	0.92	0.204	2.50	-	-	*	PeCDF	22
Total Penta-Furans	3.50e+07			30:21	0.92	94.0	2.50	-	-	*	94.2	12
Total Hexa-Furans	6.52e+07			35:37	1.34	202	2.50	-	-	*		14
Total Hepta-Furans	2.29e+07			42:23	1.06	103	2.50	-	-	*		9

Analyst: Date: 

# PROJECT REQUEST SHEET

Project #:	<u>5340</u>	Sample #:	<u>1 - 2</u>	Client Manager:	<u>BS</u>
Client:	<u>USEPA Region 5</u>		Hold Time:	<u>02/18/2009</u>	
Matrix:	<u>Soil</u>	Extraction Batch:	<u>X1652</u>	Due Date:	<u>02/26/2009</u>
Method:	<u>DLM02.0 D/F</u>		Storage:	<u>R1</u>	
SOP:	<u>SOPs: EP2C Rev.1 IP2C Rev.2</u>				

## COMMENTS/INSTRUCTIONS:

Use 1g equiv.  
Probably very high levels.

Results: 5340

5340TCDF

Extract/s located in box: Chips, Dips, Chunks, Whops

Standards: 5340

Instrument:

DB5 FAC-3

DB225 FAC-1

DB1 \_\_\_\_\_

Other \_\_\_\_\_

## Percent Solids

FAL Project: 5340

Sample ID	Chemist	Date	Wet Sample Weight (g)	Dry Sample Weight (g)	% Solids	10g Equiv
5340-001-0001-SA	GB	2/13/09	4.95	3.72	75.15	13.31
5340-002-0001-SA	↓	↓	↓	↓	↓	↓

### % Solids Summary:

#### Non-Filtered Determination

1. Place an aliquot of sample into a pre-weighed aluminum weighing boat. Use approximately two to ten grams for solid samples, approximately 10 mL for aqueous samples.
2. Record the weight.
3. Dry sample overnight at approximately 110 C.

#### Filtered Determination

1. Pre-weigh a glass fiber filter of appropriate pore size and pressure filter a sample aliquot (200-1000mL) through it.
2. Air dry the filter and record the dry weight.

#### % Solids calculation

$$\% \text{ solids} = \text{aliquot after drying}/\text{aliquot before drying} \times 100$$

- Samples containing one percent solids or less are prepared as aqueous samples.
- Samples containing greater than one percent solids prepared as solid samples.

## EXTRACTION SHEET

Project #: 5340 Extraction Date: 2009-02-18 Extraction Chemist:BM

Method/Analysis: DLM02.0 D/F

Procedure: SOX/SDS

Solvent: Toluene

Sample ID	Wet wt. (g/L)	Dry wt. (g/L)	IS		NS		CSS	
			Amt: ID:	Vial: Chemist/Witness/Date	Amt: ID:	Vial: Chemist/Witness/Date	Amt: ID:	Vial: Chemist/Witness/Date
1652-001-0001-MB	(1.00)	(1.00)	BM GN 2/18/09		NA	BM GN 2/19/09		
1652-001-0001-OPR	(1.00)	(1.00)			BM GN 2/18/09			
5340-001-0001-SA	1.33	1.00			NA			
5340-002-0001-SA	1.35	1.01		↓	↓	↓		

AX-21 Charcoal Cleaned	102708	Acetone	48233	Acid Alumina	14003MH	Hexane	48330
Hydrochloric Acid	B08505	Methanol	085079	Methylene Chloride (DCM)	48352	Silica Gel	TA1493434
Sodium Hydroxide	8120	Sodium Sulfate	48094820	Sulfuric Acid	084858	Tetradecane	072385
Toluene	48305	Water	48281	C-18 Empore Discs	320441	Cyclohexane	48149

Comments:

## CLEANUP SHEET

Project #: 5340

Method/Analysis: DLM02.0 D/F

Splits: 0 Split Date: N/A Final Volume: 20.0µL

Sample ID	Cleanup 1	Cleanup 2	Cleanup 3	RS
	<u>MS6/AA</u>	<u>CC</u>	<u>NA</u>	Amt: <u>10.0µL</u> ID: <u>080827D</u> Vial: <u>3</u> Chemist/Witness/Date
1652-001-0001-MB	<u>BM 2/19/09</u>	<u>BM 2/19/09</u>	<u>NA</u>	<u>BM GN 2/19/09</u>
1652-001-0001-OPR				
5340-001-0001-SA				
5340-002-0001-SA	<u>↓</u>	<u>↓</u>	<u>↓</u>	<u>↓</u>

Comments:

Frontier Analytical Laboratory

EPA  
1 Sample Jar

Sample Check in/Check out/Disposal Log

FAL Project No. 5340

Check in / Check out

#2 = duplicate sample, take from same jar labeled 5340-001-SA.

Sample Number(s)	Dupes	Check In	Check Out	To/From Location	Date	Time	Entire Sample Used?	Initials
1 - 2*	0	X		R1	2009-02-11	10:45:00		BM
+			X	R1	2/12/09	16:00		bb
1 (1st)			X	R1	2/13/09	16:00		bb
1 (1st)		X		R1	2/13/09	16:15		bb
1 (1st)			X	R1	2/18/09	13:10		BM
1 (1st)		X		R1	2/18/09	14:00		BM
1 (1st)			X	R1	2/25/09	15:20		GG
1 (1st)		X		R1	2/25/09	16:10		GG

pH

pH Range	Sample Number/s
pH < 4	
4 < pH < 10	NA
10 < pH	

Disposal

Sample Number/s	Over 1ppb*	Under 1ppb	From Location	Disposal Date/Time	Initials

\*Samples over 1ppb total of PCDDs/PCDFs, PCBs or PAHs require specialized disposal.  
Contact SafetyKleen to arrange disposal of any expired samples over the 1 ppb level.

# FedEx® US Airbill

Express

FedEx  
Tracking  
Number

8681 4478 7343

**1 From**

Date 2-10-09

Sender's  
Name

Company SULLIVAN INDUSTRIAL GROUP

Address 2325 VICKEROO RIVER DR.

City ELGIN

State MI

Dept/Floor/Suite/Room

Phone 517 599 1875



**2 Your Internal Billing Reference**

**3 To**

Recipient's  
Name

SAMPLE RECEIPT

Phone

Company FRONTIER ANALYTICAL LABORATORY

Recipient's  
Address

5172 HILLSONG CIRCLE

Dept/Floor/Suite/Room


Address

To request a package be held at a specific FedEx location, print FedEx address here.

City EL DORADO HILLS

State CA

ZIP 95762



8681 4478 7343

Form  
ID No.

0200

Recipient's Copy

**4a Express Package Service**

FedEx Priority Overnight  
Next business morning\* Friday  
shipments will be delivered on Monday  
unless SATURDAY Delivery is selected.

FedEx 2Day  
Second business day\* Thursday  
shipments will be delivered on Monday  
unless SATURDAY Delivery is selected.

FedEx Envelope rate not available. Minimum charge: One-pound rate.

FedEx Standard Overnight  
Next business afternoon\* Saturday  
Delivery NOT available.

FedEx Express Saver  
Third business day\* Saturday  
Delivery NOT available.

**Packages up to 150 lbs.**

FedEx First Overnight  
Earliest next business morning  
delivery to select locations.\*  
Saturday Delivery NOT available.

\* To most locations.

**4b Express Freight Service**

FedEx 1Day Freight\*  
Next business day\*\* Friday  
shipments will be delivered on Monday  
unless SATURDAY Delivery is selected.

FedEx 2Day Freight  
Second business day\* Thursday  
shipments will be delivered on Monday  
unless SATURDAY Delivery is selected.

**Packages over 150 lbs.**

FedEx 3Day Freight  
Third business day\*\* Saturday  
Delivery NOT available.

\*\* To most locations.

\* Call for Confirmation:

**5 Packaging**

FedEx  
Envelope\*

FedEx Pak\*  
Includes FedEx Small Pak,  
FedEx Large Pak, and FedEx Sturdy Pak.

FedEx  
Box

FedEx  
Tube

Other  
\* Declared value limit \$500.

Include FedEx address in Section 3.

SATURDAY Delivery  
Not available for  
FedEx Standard Overnight,  
FedEx First Overnight, FedEx Express  
Saver, or FedEx 3Day Freight.

HOLD Weekly  
at FedEx Location  
Not available for  
FedEx First Overnight.

HOLD Saturday  
at FedEx Location  
Available ONLY for FedEx Priority  
Overnight and FedEx 2Day  
to select locations.

**Does this shipment contain dangerous goods?**

No  Yes  
As per attached  
Shipper's Declaration.

Yes  
Shipper's Declaration  
not required.

Dry Ice  
Dry Ice, 9 UN 1845 x kg  
 Cargo Aircraft Only

**7 Payment Bill to:**

Sender  
Acct. No. in  
Section 1 will  
be billed.

Recipient

Third Party

Credit Card

Cash/Check  
Obtain Recip.  
Acct. No.

Total Packages Total Weight Total Declared Value\*

\$ .00

\*Our liability is limited to \$100 unless you declare a higher value. See back for details.

**8 Residential Delivery Signature Options**

No Signature  
Required  
Package may be left  
without obtaining a  
signature for delivery.

Direct Signature  
Someone at recipient's  
address may sign for  
delivery. *Fee applies.*

Indirect Signature  
If no one is available at  
recipient's address, someone  
at a neighboring address may  
sign for delivery. *Fee applies.*

520

Rev. Date 10/06 Part #158281 ©1994-2006 FedEx PRINTED IN U.S.A. SRY

Credit Card Auth.

L

000275 of 000279



**USEPA Contract Laboratory Program  
Organic Traffic Report & Chain of Custody Record**

Analytic Number: P0# EP08W001509

1. Case No.: <u>PO# EP08W001509</u>	2. Date Shipped: <u>2/11/2009</u>	3. Chain of Custody Record	Sampler Signature: <u>Jerry Beach</u>	4. For Lab Use Only						
DAS No.: <u>SDG</u>	Carrier Name: <u>FCD BX</u>	Relinquished By: <u>Jerry</u>	Received By: <u>Jerry Beach</u>	Lab Contract No: <u>EP08W001509</u>						
SDG No.: <u>E3YCA</u>	Airbill: <u>8432 05044404</u>	(Date/Time) <u>2/11/09 1600</u>	Date/Time) <u>2/11/09 1600</u>	Unit Price: <u>\$100</u>						
Shipped To: <u>Frontier Analytical Lab</u> <u>5172 Hillsdale Circle</u> <u>El Cajon, CA 92112</u>		Transfer To:								
5. ORGANIC SAMPLE No.	6. MATRIX/ SAMPLER	7. TYPE	8. ANALYSIS/ TURNAROUND	9. TAG No./ PRESERVATIVE/ Bottles	10. STATION LOCATION	11. SAMPLE COLLECT DATE/TIME	12. INORGANIC SAMPLE No.	13. FOR LAB USE ONLY Sample Condition On Receipt		
E3YCA	SOLV PRODUCT	GAB	Dixie/ Farrow	1 80707JUN	January 13 2009	13 30	5340-001-SA First	good, 2/11/09 10:45		
E3YCA	SOLV PRODUCT	GAB	Dixie/ Farrow	1 80707JUN	1/13/2009	13 30	5340-002-SA Last	good, 2/11/09 10:45		
14. Shipment for Case Complete? <input checked="" type="checkbox"/> 15. Sample(s) to be used for laboratory QC: <input checked="" type="checkbox"/> 16. Additional Sampler Signature(s): <input checked="" type="checkbox"/>									17. Cooler Temperature Upon Receipt: <u>18°C</u>	18. Chain of Custody Seal Number: <u>18C</u>
19. Analysis Key: <input checked="" type="checkbox"/> Type: Comp, Grab (from Box 7)									20. Custody Seal Intact? <input checked="" type="checkbox"/>	21. Shipment Iced? <input checked="" type="checkbox"/>

22. TR Number:  
00276

23. R provides preliminary results. Requests for preliminary results will increase analytical costs.  
Send Copy to: Sample Management Office, Attn: Heather Bauer, 15000 Conference Center Dr., Chantilly, VA 20151 Phone 703/818-4200 Fax 703/995-4415

**LABORATORY COPY**

Page 1 of 6

CB CONGENER  
SAMPLE LOG-IN SHEET (DC-1)

Lab Name <u>Frontier Analytical Laboratory</u>		Page <u>1</u> of <u>1</u>			
Received By (Print Name) <u>TOM CRABTREE</u>	Log-in Date <u>2/11/09</u>				
Received By (Signature) <u>TOM CRABTREE</u>					
Contract No. <u>EPA08W001564</u>	TO No.				
Case No. <u>EPA05-2008-01</u>	Sample Delivery Group No.				
Remarks:	EPA Sample #	Corresponding		Remarks: Condition of Sample Shipment, etc.	
	Sample Tag #	Assigned Lab #			
	<u>E3YLC9</u>		<u>5340-001-SA</u>		<u>Good, First</u>
	<u>E3YD8</u>		<u>5340-002-SA</u>		<u>Good, Last</u>
Sample Transfer					
Fraction <u>#1</u>	Fraction <u>#1</u>				
Area # <u>SC-frap</u>	Area # <u>Pop-Int</u>				
By <u>BM</u>	By <u>BM</u>				
On <u>2/18/09</u>	On <u>2/19/09</u>				

\* Contact SMO and attach record of resolution.

Reviewed By <u>JJS</u>	Logbook No.
Date <u>3/5/09</u>	Logbook Page No.

**Gabby Navarro**

---

**From:** Ben Miller [benm@frontieranalytical.com]  
**Sent:** Thursday, February 12, 2009 2:41 PM  
**To:** 'vandonsel.terese@epamail.epa.gov'  
**Cc:** 'Kathy Zipp'; 'brads@frontieranalytical.com'  
**Subject:** 5340 Sample Receipt EP08W001564 Due: 2/26/2009

Dear Ms. Donsel,

This letter confirms the receipt of two soil samples at Frontier Analytical Laboratory on February 11, 2009. We received two soil samples and a total of 1 bottle. The bottle arrived in good condition. The samples will be extracted and analyzed using EPA Method DLM02.0 for tetra through octa chlorinated dibenzo dioxins and furans. These samples are associated with EPA Order/Contract No.: **EP08W001564**. This project has a turnaround time of fourteen calendar days at \$900/sample. The following e-mail is a PDF copy of the completed Chain-of-Custody/Traffic Report. The samples have been assigned a Frontier Analytical Laboratory project number **5340** and SDG Number **SDG E3YC9** In order to expedite inquiries, please reference our project number.

These samples represent the first samples in SDG Number **SDG E3YC9**.

**Please note that the samples were received at 18 degrees C which is outside the method recommended temperature range of 0 – 4 degrees C. We will proceed with analysis unless instructed otherwise.**

If you have any questions regarding this project, please feel free to contact us at your convenience.

Thank you for choosing Frontier Analytical Laboratory for your analytical testing needs.

Sincerely,

Benjamin Miller

Sample Control  
Frontier Analytical Laboratory  
5172 Hillsdale Circle  
El Dorado Hills, CA 95762  
Tel: (916) 934-0900  
Fax: (916) 934-0999  
[www.frontieranalytical.com](http://www.frontieranalytical.com)

*The materials transmitted by this electronic mail are confidential, are only for the use of the intended recipient, and may also be subject to applicable privileges. Any dissemination, distribution, or copying of this communication is strictly prohibited. If you have received this communication in error, please immediately notify the sender. Please also remove this message from your hard drive, diskette, and any other storage device.*

**Gabby Navarro**

---

**From:** Sample Control [samplecontrol@frontieranalytical.com]  
**Sent:** Thursday, February 12, 2009 2:44 PM  
**To:** 'vandonsel.terese@epamail.epa.gov'  
**Cc:** 'brads@frontieranalytical.com'; 'Kathy Zipp'; 'gabbyn@frontieranalytical.com'  
**Subject:** Dioxin Contract No. EP08W001564

Dear Ms. Donsel,

We received one 8oz jar for sample E3YC9 and E3YDO in good condition on February 11, 2009. We found the following discrepancies:

- The Chain of Custody has hand written Case No. of PO# EP08W001564; however our records show a Case No. of EPA05-2008-01. Is the Case No. the same as the DAS No.?
- We didn't receive sample tags with the samples. Does your region require sample tags?
- Does your region require a Delivery of Analytical Services (DAS) Sample Delivery Group (SDG) Sheet and if so can you send us EPA Region 5's form?
- Does your region require a TR Number? We didn't receive a number with the Chain of Custody (COC).
- Please note that the Airbill number on the COC of 863205044404 does not match the Airbill number received of 868144787343.

Please feel free to contact us at your convenience.

Thank you for choosing Frontier Analytical Laboratory for your analytical testing needs.

Sincerely,  
Gabby Navarro

Sample Control  
Frontier Analytical Laboratory  
5172 Hillsdale Circle  
El Dorado Hills, CA 95762  
Tel: (916) 934-0900  
Fax: (916) 934-0999  
[www.frontieranalytical.com](http://www.frontieranalytical.com)

*The materials transmitted by this electronic mail are confidential, are only for the use of the intended recipient, and may also be subject to applicable privileges. Any dissemination, distribution, or copying of this communication is strictly prohibited. If you have received this communication in error, please immediately notify the sender. Please also remove this message from your hard drive, diskette, and any other storage device.*