## **APPENDIX C**



# **QUALITY ASSURANCE PROJECT PLAN**

**FOR** 

### **REMOVAL ACTIVITIES**

AT THE

#### TOLEDO TIE TREATMENT SITE

LOCATED AT

ARCO INDUSTRIAL PARK TOLEDO, OHIO

> FEBRUARY 1998 (Revised April 1998)

> > Prepared Fors

KERR-McGEE CHEMICAL, LLC KERR-McGEE CENTER OKLAHOMA CITY, OKLAHOMA 73125

Prepared By:

HULL & ASSOCIATES, INC. 2726 MONROE STREET TOLEDO, OHIO 43606

### TABLE OF CONTENTS

•	, in the second of the second	Page
C t A	INTRODUCTION	1
C.1.0	NTRODUCTION	1
	C.1.1 General	1
	C.1.2 Objective	
C.2.0	PROJECT ORGANIZATION AND RESPONSIBILITIES	3
	C.2.1 General	3
	C.2.2 Field Operations Personnel	3
	C.2.2.1 Project Manager (PM)	3
•	C.2.2.2 Quality Assurance Officer (QAO)	
:(	C.2.2.3 Field Operations Coordinator (FOC)	4
	C.2.3 Laboratory Personnel	Δ
	C.2.3.1 Laboratory Director	Δ
	C.2.3.2 Laboratory Project Director	5
	C.2.3.3 Laboratory Quality Assurance Officer	
C.3.0	QA OBJECTIVES FOR MEASUREMENT	6
	·	
	C.3.1 General	
	C.3.2 Precision	
	C.3.3 Accuracy	6
	C.3.4 Representativeness	
	C.3.5 Completeness	
	C.3.6 Comparability	7
C 4 A		O
C.4.0	SAMPLING	ð
	C.4.1 General	Q
	C.4.2 General Sampling Procedures	
	C.4.2 General Sampling Procedures	0
	C.4.2.1 Sample Containers and Preservatives	8
	C.4.2.2 Sample Labeling	
	C.4.2.3 Sample Identification Numbers	
	C.4.2.4 Sampling Equipment Preparation and Decontamination	
•	C.4.2.5 Sample Storage and Transportation	
	C.4.2.6 Field Notes	

### **TABLE OF CONTENTS (cont.)**

	, · · · · ·	Page
	C.4.2.7 Chain-of-Custody	11
	C.4.2.8 Field Sampling Equipment List	11
	C.4.2.9 Sampling Quality Control	
	C.4.3 Site-Specific Sampling Procedures	12
C.5.0	SAMPLE CUSTODY	13
	C.5.1 General	13
	C.5.2 Chain-of-Custody.	
	C.5.3 Laboratory Custody Procedures	
	C.5.4 Laboratory Documentation.	
C.6.0	CALIBRATION PROCEDURES AND FREQUENCY	16
	C.6.1 General	16
	C.6.2 Field Equipment Calibration Procedures	
	C.6.3 Laboratory Instrumentation Calibration Procedures	
C.7.0	ANALYTICAL PROCEDURES	18
C.8.0	DATA REDUCTION, REVIEW, REPORTING AND VALIDATION	19
	C.8.1 Data Reduction	19
	C.8.2 Data Review	
	C.8.3 Data Reporting.	20
	C.8.4 Data Validation	20
C.9.0	INTERNAL QUALITY CONTROL	21
	C.9.1 General	. 21
	C.9.2 Blank Samples	21
	C 0 2 1 Mathed Planks	21
	C.9.2.1 Method Blanks	
	C.9.2.2 Container Blank	
	C.9.2.3 Holding Blanks	22
	C.9.3 Reference Standards	22
	C.9.4 Analytical Spikes	23

### TABLE OF CONTENTS (cont.)

		Page
	C.9.4.1 Matrix Spike	23
	C.9.4.2 Surrogate Spike	
	C.9.5 Replicate Analysis.	23
	C.9.6 Calibration Check Standards	
	C.9.7 Internal Standards	
C.10.0	D PERFORMANCE AND SYSTEM AUDITS	25
C.10.		•
	C.10.1 General	25
	C.10.2 Performance Audits	25
	C.10.3 System Audits	25
	C.10.4 Field Audits	26
C.11.0	PREVENTIVE MAINTENANCE	27
		,
C.12.0	DATA PRECISION, ACCURACY, AND COMPLETENESS	28
	C.12.1 Precision	28
	C.12.2 Accuracy	28
· .	C.12.3 Completeness	
C 13 (	CORRECTIVE ACTION	. 20
C.13.0		
	C.13.1 General	30
	C.13.2 Laboratory Corrective Action	
		20
	C.13.2.1 Bench Level	
	C.13.2.2 Management Level	
	C.13.2.3 Receiving Level	31
•	C.13.3 Field Corrective Action	31
C 14 (1	QUALITY ASSURANCE REPORTS TO MANAGEMENT	32
<b>∵.</b> ≀7.⊍	QUALITY ADDURANCE REFORM ONTO TO MANAGEMENT	2
•	C.14.1 Internal Reporting	32
	C.14.1 Internal Reporting C.14.2 Additional Reporting	32
C.15.0	REFERENCES	33

### **TABLE OF CONTENTS (cont.)**

#### LIST OF TABLES

Parameters List
Sampling Containers and Preservatives
Sample Types and Codes
Field Equipment List
Analytical Procedures

### LIST OF FIGURES

Figure C1 Project Organization Chart

#### LIST OF ATTACHMENTS

Attachment A	Lancaster Laboratories' Quality Assurance Plan
Attachment B	Field Data Sheets

TOLEDO TIE TREATMENT SITE QUALITY ASSURANCE PROJECT PLAN FINAL

PAGE 1 of 33

C.1.0 INTRODUCTION

C.1.1 General

Kerr-McGee Chemical, LLC (Kerr-McGee) was issued a Unilateral/Administrative Order (UAO) dated

December 24, 1997, pursuant to Section 106(a) of the Comprehensive Environmental Response,

Compensation, and Liability Action (CERCLA) pertaining to the Toledo Tie Treatment Site located in

Toledo, Ohio (Site). Preparation of this document is in accordance, to the maximum extent practicable,

with the provisions of the UAO, Docket No. V-W-'98-C-444 which requires Kerr McGee to prepare

work plans necessary to conduct a time-critical removal action. A more detailed description of the Site

and the proposed work activities is included in Section 1.0 of the Work Plan (HAI Document No.

PWM001D.002).

This Quality Assurance Project Plan (QAPP) has been prepared for personnel representing Kerr-McGee,

the U.S. Environmental Protection Agency (EPA), and HAI field personnel conducting the investigation.

The QAPP is intended to provide the quality assurance and quality control guidelines for the activities

described in the Work Plan and the Field Sampling and Analysis Plan (FSAP) which is Appendix A of

the Work Plan.

The analytical subcontractor for this project will be Lancaster Laboratories. All analyses will be

performed at their Lancaster, Pennsylvania facility. Lancaster's QAPP has been included as Attachment

A of this document. If an alternate laboratory is required, its QAPP will be submitted as an addendum

to this plan.

C.1.2 Objective

The objective of this plan is to document the procedures that will be used to collect sufficient data of

known quality for the time-critical Removal Action. This plan is intended to be used as supplemental

guidance to the Work Plan and Field Sampling and Analysis Plan. The quality assurance measures for

HULL & ASSOCIATES, INC. TOLEDO, OHIO

TOLEDO TIE TREATMENT SITE QUALITY ASSURANCE PROJECT PLAN FINAL PAGE 2 of 33

the analytical program will be in accordance with the appropriate U.S. EPA methods, good laboratory practices, and the laboratory's quality assurance program (located in Attachment A).

Table C.1 provides the anticipated list of the parameters to be analyzed and the accompanying methods. Reporting limits for all parameters will be in accordance with the provisions of the analytical methods used and good laboratory practices. Reporting limits may vary between samples as they can be affected by sample matrix, dilutions, and other interferences. If additional analytical methods are required, this decument will be revised accordingly by addendum.

#### C.2.0 PROJECT ORGANIZATION AND RESPONSIBILITIES

#### C.2.1 General

The organizational chart for the project is presented in Figure 1. The individual who is primarily responsible for Quality Assurance (QA) will be the Quality Assurance Officer (QAO), who reports directly to the Project Manager. The responsibilities of key individuals are outlined below.

#### C.2.2 Field Operations Personnel

#### C.2.2.1 Project Manager (PM)

The PM for this project will be Scott Lockhart, P.E. of HAI. Mr. Lockhart will be responsible for the overall development and management of the project and is the communication link between project personnel, and any applicable regulatory agencies. Duties and responsibilities of the PM will be to:

- 1. administrate and supervise all phases of the project;
- 2. ensure project objectives are met within financial and time constraints;
- 3. work with the QAO and field personnel to plan and conduct project operations, progress meetings, etc.; and
- 4. review progress reports and analytical reports prior to being issued.

#### C.2.2.2 Quality Assurance Officer (QAO)

The QAO for this project will be Kevin Wildman of HAI. Mr. Wildman will be responsible for the adherence to the QAPP. Duties and responsibilities of the QAO will be to:

- 1. establish QA/QC procedures required for the project;
- 2. evaluate data quality and maintain QC records;
- 3. provide a communication link between project personnel and the laboratory;

- 4. monitor the progress of the field sampling personnel and enforce provisions of this plan; and
- 5. stop work at any time that the QAPP is not being adhered to, or if the quality of the results are jeopardized by the work in progress; once work is stopped, only the PM can restart activities.

#### C.2.2.3 Field Operations Coordinator (FOC)

The FOC for this project will be Mr. Jeff Arp. Mr. Arp will be responsible for overseeing the day-to-day conduct of project activities. Duties and responsibilities of the FOC will be to:

- 1. ensure the sampling activities are conducted in a manner that follows the procedures outlined in this plan and the Work Plan;
- 2. coordinate the sampling activities with the PM, QAO, and field personnel;
- 3. oversee the use, maintenance and operation of sampling equipment; and
- 4. report daily activities, problems, etc. to QAO and PM.

#### C.2.3 Laboratory Personnel

The laboratory will have its own project organization with responsibilities similar to that of the field operations personnel.

#### C.2.3.1 Laboratory Director

The Laboratory Director will be primarily responsible for the overall operation of the laboratory including all samples analyzed and data reported. The Laboratory Director will also be responsible for initiating corrective action measures when analytical data do not meet the requirements of this plan or the laboratory's QAPP.

TOLEDO TIE TREATMENT SITE QUALITY ASSURANCE PROJECT PLAN FINAL PAGE 5 of 33

#### C.2.3.2 Laboratory Project Director

The Laboratory Project Director will be the primary communication link between the laboratory and HAI's QAO. The Laboratory Project Director will be responsible for relating any special needs of the field operations personnel to the laboratory. The Laboratory Project Director will also provide the final review of all data packages before reporting results.

#### C.2.3.3 Laboratory Quality Assurance Officer

The Laboratory QAO will be primarily responsible for implementing the laboratory's QAPP within the laboratory and monitoring compliance with the laboratory's QAPP. The Laboratory QAO's duties will also include: conducting audits, reviewing all QC sample recoveries, reporting problems to Laboratory Director for corrective action, and other laboratory-related activities.

HULL & ASSOCIATES, INC. TOLEDO, OHIO

FINAL PAGE 6 of 33

C.3.0 QA OBJECTIVES FOR MEASUREMENT

C.3.1 General

Data quality objectives for measurement during this project will be addressed in terms of precision,

accuracy, representativeness, completeness, and comparability (PARCC parameters). The collection of

data used in this project will require that the sampling and analysis be performed using standard methods,

with properly operated and calibrated equipment, and conducted by trained personnel.

C.3.2 Precision

Precision is the determination of the reproducibility of measurements under a given set of conditions, or

a quantitative measure of the variability of a group of measurements compared to their average value.

Precision of analytical results will be based upon laboratory replicate analyses. Precision is reported as

Relative Percent Difference (%RPD). Precision goals for the parameters to be analyzed will be in

accordance with the provisions of the U.S.EPA methods used for analysis.

C.3.3 Accuracy

Accuracy is defined as the degree of agreement of a measurement or average of measurements with an

accepted reference or true value. Sampling accuracy can be assessed by evaluating the results of field

and trip blanks. Analytical accuracy is assessed by percent recovery of analytical spikes and reference

standards. Accuracy goals for parameters to be analyzed will be in accordance with the provisions of

the U.S.EPA methods used for analysis.

C.3.4 Representativeness

Representativeness expresses the degree to which sample data accurately and precisely represent a

characteristic of a population, parameter variations at a sampling point, and/or environmental condition.

Representativeness is best achieved by insuring that sampling locations are properly selected and

sufficient number of samples are collected.

HULL & ASSOCIATES, INC. TOLEDO, OHIO

TOLEDO TIE TREATMENT SITE QUALITY ASSURANCE PROJECT PLAN

FINAL PAGE 7 of 33

Functional data are provisional on meeting the criteria for representativeness. The QA goal will be to

have all samples and measurements representative of the media sampled; representative samples are

contingent on proper selection of sampling techniques, location, and number of samples collected to

represent the media. Also, the aliquots taken for analysis should be representative of the sample

received.

C.3.5 Completeness

Completeness expresses the measure of confidence with which the data, resulting from a data collection

activity, meets the specific objectives of the activity. It is the amount of valid data obtained from a

measurement system compared to the amount that was expected to be obtained under correct normal

operating conditions. While efforts will be made to have all generated data be valid data (complete),

some samples may be lost or broken in transit or at the laboratory. In addition, the results may not be

acceptable based on laboratory QC requirements. It will be the goal of this project to have 95% of the

data generated to be complete data.

C.3.6 Comparability

Comparability expresses the measure of confidence in which data sets can be considered equivalent with

regard to the measurement of a specific parameter and/or groups of parameters. The ability to compare

data sets is particularly critical when a set of data for a specific parameter is compared to historical data

for determining trends. By using standard sampling techniques, analytical methodologies, and reporting

units, the comparability of data will be ensured.

#### C.4.0 SAMPLING

#### C.4.1 General

The purpose of this section is to detail the general sampling procedures that will be used to collect the data required to complete this project. The sampling efforts shall be uniform and follow specific protocols to be considered relevant to the project. Additional information is included in the Work Plan.

#### **C.4.2** General Sampling Procedures

#### C.4.2.1 Sample Containers and Preservatives

Sample containers will consist of I-Chem 200 series (or equivalent) glass or plastic bottles and will be provided and prepared by the laboratory prior to sampling efforts. The laboratory will also provide any required preservatives. Table C.2 lists the containers, preservatives, and holding times for parameters that will be analyzed during this project.

#### C.4.2.2 Sample Labeling

All sample containers will be labeled at the time of sampling. Each label will be completed with the required information and then secured to the container with transparent packing tape to prevent accidental loss or damage. Required information on the sample label includes: project number, sample identification number, date, time, analyses, and sampler's initials. Also, any preservatives or special handling instructions will be clearly displayed on the label.

#### **C.4.2.3 Sample Identification Numbers**

All samples collected will be issued a unique Sample Identification Number (SIN) to aid in tracking and record keeping. SINs will be assigned by the QAO, and given to the FOC prior to sampling efforts. SINs will consist of four parts separated by hyphens. The first part of a SIN will be HAI's project number which is "PWM001."

The second part of a SIN represents the type of sample collected. Valid sample types and codes for this project are listed in Table C.3.

HULL & ASSOCIATES, INC.

APRIL 1998
TOLEDO, CHIO
PWM001D.001

**APRIL 1998** 

The third part of a SIN indicates the frequency of sampling, which sample of a similar sample type is referenced, or the interval from which a soil boring sample was collected. Numbering of the above sample types shall be as follows:

- 1. The majority of investigative samples (e.g., sediment, surface soil, surface water, etc.) will use a three digit number to indicate that a new sample location is being referenced. Numbers will start at 001 and increase by one for every new sample of that sample type.
- 2. Monitoring wells will use a three digit number to indicate the number of times the well has been sampled (i.e., the first time a monitoring well is sampled, the number will be 001, and the second time the number will be 002). When collecting a duplicate sample. the sample frequency number will be the same as the original sample; however, an "A" will be added to identify this as a duplicate sample
- 3. Soil boring samples are numbered slightly different than described above. The second segment of a soil boring sample number will designate the soil boring number, while the third segment will indicate the depth interval from which the sample was obtained. The first sample collected from a boring will be designated SS1, and the second SS2, etc. The corresponding depths of each soil sample will be clearly identified on the Soil Boring Log.
- 4. Field blanks and trip blanks will use a six digit number to indicate the date that the sample was collected. Format of this numbering system will be month, date, and year. An example for a blank collected on May 1, 1998 would be 050198.

The fourth and final part of a SIN will be a four digit code to identify the person responsible for collecting the sample. This code will consist of the first letter of the office from which the sampler originates followed by the individuals employee number (i.e., D153 would be employee number 153 from the Dublin office).

Designated office codes are:

D - HAI Dublin, OH T - HAI Toledo, OH M - HAI Mason, OH W - HAI Warrensville Heights, OH

PAGE 10 of 33

The SIN system described above is very important in the tracking and record keeping of the large number

of samples to be collected for this project. For this reason, the SIN system will not be deviated from

without authorization of the QAO. Any questions regarding the SIN system will be directed towards

the QAO.

C.4.2.4 Sampling Equipment Preparation and Decontamination

Sampling equipment to be reused will be thoroughly decontaminated between sampling locations and at

the beginning and end of each day. To decontaminate the equipment it will be washed with a mild non-

phosphatic soap and thoroughly rinsed with distilled water. HAI Standard Operating Procedure (SOP)

No. F1000 (refer to Attachment A of the FSAP) provides a more detailed description of decontamination

procedures. If complete cleaning of any piece of sampling equipment is not possible, then it will be

discarded and replaced with a clean article.

C.4.2.5 Sample Storage and Transportation

Field samples will be placed in portable coolers on ice immediately following sample collection and

remain on ice until being delivered to the laboratory. Ice will be double bagged to prevent leakage and

possible water damage to samples, sample labels, and documentation. Any samples not placed on ice

immediately upon collection will be discarded, and a new sample will be collected.

C.4.2.6 Field Notes

General field notes will be recorded in waterproof surveyors notebooks using indelible ink. In addition

to the field notebooks, certain activities will require the completion of data sheets. A Soil Boring Log

(see Attachment B) will be completed for each soil boring/monitoring well installed. A Groundwater

Monitoring Well Sampling Data Sheet (see Attachment B) must be completed for each monitoring well

sampled. When weather conditions prohibit the completion of data sheets in the field, data may be

recorded in field notebooks and then transferred to data sheets at the end of the day.

Additionally, a Daily Field Report (Attachment B) will be completed at the end of the day summarizing

the day's activities and observations. Copies of the documentation will be forwarded to HAI's Dublin,

TOLEDO TIE TREATMENT SITE QUALITY ASSURANCE PROJECT PLAN

PAGE 11 of 33

Ohio office weekly. If copies of previous work are required, then arrangements will be made with the

QAO.

Field notebooks, field data sheets, or daily field reports will not be obscured, destroyed, or discarded,

even if it contains errors or is illegible. Corrections will be made by drawing a single line through the

error and writing in the correct information. Corrections will be dated and initialed by the person making

the correction.

C.4.2.7 Chain-of-Custody

The chain-of-custody is discussed in Section C.5.0 of this plan.

C.4.2.8 Field Sampling Equipment List

Table C.4 is a list of the general field sampling equipment that will be available on-site. The field

analysis equipment will be calibrated in accordance with the manufacturer's recommendations and this

plan.

C.4.2.9 Sampling Quality Control

Several sampling quality control measures will be necessary to assess the integrity of samples collected.

These measures include the use of field blanks and trip blanks to locate possible sources of sample

contamination.

The number of field blanks (e.g., equipment/rinseate blanks) analyzed for a class of compounds will be

"equal to" ten percent of the total samples analyzed, for that class, with a minimum of one per day. Field

blanks will be collected by running laboratory prepared deionized water through sample collection

equipment and preserved according to Table 2. Field blanks will be analyzed for the same parameters

as the field samples. It is the samplers responsibility to collect the appropriate number of field blanks

for the day's sampling efforts.

TOLEDO TIE TREATMENT SITE QUALITY ASSURANCE PROJECT PLAN FINAL

PAGE 12 of 33

One trip blank per shipping container (e.g., cooler) will be required. Trip blanks are only necessary for samples requiring volatile organic analyses. Trip blanks will be prepared in the laboratory, prior to sampling efforts, using laboratory-prepared deionized water and preserved using the same procedures as the samples. Trip blanks must accompany sample containers during sample collection and transportation. When sampling groundwater and surface water, a field duplicate sample will be collected. The minimum frequency of field duplicate sample collection is one per every ten investigative samples. It is required that a field duplicate be collected every day during groundwater sampling events. A new

C.4.3 Site-Specific Sampling Procedures

Site-specific sampling procedures are presented in the Field Sampling and Analysis Plan.

field duplicate will be required if the members of the sampling team change during the day.

#### **C.5.0 SAMPLE CUSTODY**

#### C.5.1 General

The intention of chain-of-custody (COC) procedures is to document in a legally defensible manner the transfer of custody of each sample from collection through analysis. Additional information regarding COC procedures is presented in HAI SOP No. F3014 (refer to Attachment A of the FSAP).

#### C.5.2 Chain-of-Custody

The importance of COC cannot be overstated. This documentation records the history of the samples' custody from acquisition to ultimate disposal. Samples collected may be used as legal evidence. As such, the hand-to-hand custody from the point of collection to delivery at the laboratory must be clearly documented. The National Enforcement Investigations Center (NEIC) of the U.S. EPA defines custody as:

- 1. the sample is in your physical possession;
- 2. the sample is within view after being in your physical possession;
- 3. the sample was in your possession and then you locked or sealed it to prevent tampering; and/or
- 4. the sample is placed in a designated secure place with limited access to authorized personnel only.

A COC form (see Attachment B) must accompany every shipping container. Each COC form will be filled out in triplicate. Information required on the COC form includes:

- 1. client information;
- 2. project information;
- samplers' names;
- sample identification numbers;
- 5. date and time of collection;
- 6. type of sample (grab or composite);

TOLEDO TIE TREATMENT SITE QUALITY ASSURANCE PROJECT PLAN FINAL

PAGE 14 of 33

7. matrix or matrices;

8. sample description;

9. number of containers:

10. requested analyses;

11. remarks (preservatives); and

12. signatures of anyone relinquishing or accepting custody.

Field samplers will be responsible for the care and custody of the samples collected until the samples are

transferred or dispatched properly. After completing a sampling event, sample custody will be

transferred to a designated person who will maintain custody of samples until they are dispatched to the

laboratory.

If samples are to be delivered to laboratory via a courier, a COC will be signed over to the courier. The

courier will keep COC forms until relinquishing custody to the laboratory. One copy of the COC form

will be retained before there is a transfer of custody to the courier. Evidence tape or custody seals will

be placed so that when the coolers are opened the seals will be broken. Transparent tape will be used

to guarantee that the seals are not accidentally removed or destroyed.

If samples will be delivered to the laboratory via commercial carrier, then the COC forms will be placed

in a watertight, Ziploc bag and taped to the inside lid of the sample cooler. Evidence tape or custody

seals will be placed so that when the coolers are opened the seals will be broken, transparent tape will

be used to guarantee that the seals are not accidentally removed or destroyed.

C.5.3 Laboratory Custody Procedures

Samples will be received in an area designated for sample receipt and storage. Upon receipt, each

sample will be assigned a unique laboratory sample identification number. This number, along with the

date received and general description, will be recorded in the laboratory's master log. HAI's QAO will

be immediately notified if there are any problems with the samples received (e.g., discrepancies between

HULL & ASSOCIATES, INC TOLEDO, OHIO

TOLEDO TIE TREATMENT SITE QUALITY ASSURANCE PROJECT PLAN

PAGE 15 of 33

COC and samples submitted, breakage, improper preservation, etc.). Additional information regarding

laboratory custody procedures is presented in Attachment A.

**C.5.4** Laboratory Documentation

Workbooks, bench sheets, instrument logbooks, and instrument printouts are used to trace the history of

samples through the analytical process and to document and relate important aspects of the work,

including the associated quality controls. All logbooks, bench sheets, instrument logs, and instrument

printouts are part of the permanent record of the laboratory. Laboratory notebooks will be periodically

reviewed by the Laboratory Section Heads for accuracy, completeness, and compliance with the

Laboratory Quality Assurance Program Plan. Completed workbooks and instrument logbooks will be

submitted to the Laboratory Director for storage.

The laboratory's documentation procedures are presented in Attachment A. In general, good laboratory

practices require that the following (or equivalent) procedures be used. Each page, or as required, each

entry will be dated and initialed by the analyst when the record is made. Errors in entry will be crossed

out in indelible ink with a single stroke. The use of white-out, obliterating, or writing directly over the

erroneous entry will be prohibited. All corrections will be initialed by the individual making the

correction.

HULL & ASSOCIATES, INC. TOLEDO, OHIO

C.6.0 CALIBRATION PROCEDURES AND FREQUENCY

C.6.1 General

This section details the calibration procedures and frequency for both the field and laboratory

instrumentation that will be used during this project. Materials used for instrument calibration will be

obtained through the U.S. EPA Pesticide and Industrial Chemicals Repository, or a suitable commercial

source.

**C.6.2** Field Equipment Calibration Procedures

Equipment to be used during the field sampling will be examined to certify that it is in operating

condition. This includes checking the manufacturer's operating manual and instructions for each

instrument to ensure that all maintenance requirements are being observed.

Instruments and equipment used to gather, generate, or measure environmental data will be calibrated

with sufficient frequency and in such a manner that accuracy and reproducibility of results are consistent

with the manufacturer's specifications. The calibration of field instruments will be in accordance with

the manufacturer's specifications. Frequency of calibration will be dictated by field conditions,

instrument response, and the manufacturer's specifications. At a minimum all instruments will be

calibrated at the beginning of each day and after any extended breaks (e.g. lunch).

C.6.3 Laboratory Instrumentation Calibration Procedures

Calibration of laboratory equipment will be based on approved written procedures. Records of

calibration, repairs, or replacement will be filed and maintained by the laboratory's QAO. These records

will be filed at the location where the work is performed and will be subject to QA audits. For all

instruments, the laboratory will maintain a factory-trained repair staff with in-house spare parts or will

maintain service contracts with vendors.

TOLEDO TIE TREATMENT SITE QUALITY ASSURANCE PROJECT PLAN FINAL PAGE 17 of 33

The laboratory will operate in strict accordance with the requirements of the U.S.EPA methods used. Any proposed deviations from calibration procedures will be submitted by the laboratory for approval in the form of a SOP.

TOLEDO TIE TREATMENT SITE QUALITY ASSURANCE PROJECT PLAN FINAL PAGE 18 of 33

#### C.7.0 ANALYTICAL PROCEDURES

The analytical methods which will be employed in this project are summarized in Table C.5. All analytical procedures will be conducted in accordance with the specified U.S.EPA methods. In the event that additional procedures are deemed necessary, the appropriate modifications will be made to this QAPP by revision or addendum.

TOLEDO TIE TREATMENT SITE QUALITY ASSURANCE PROJECT PLAN

FINAL.

PAGE 19 of 33

C.8.0 DATA REDUCTION, REVIEW, REPORTING AND VALIDATION

C.8.1 Data Reduction

Analytical results will be reduced to the concentration units using the equations specified in the analytical

procedure. Appropriate blank corrections will be applied in all cases. Calculations will be checked by

senior laboratory staff.

C.8.2 Data Review

Each laboratory section will provide extensive data review prior to reporting results. In general, there

are three levels of data review.

The analyst will be responsible for primary review of data generated from sample analysis. If recoveries

of all quality control samples are within the method specified tolerances then the data will be presented

to data review groups for secondary review. If recoveries of any quality control samples exceed

specified tolerances, affected samples will be re-analyzed.

Secondary review will be conducted by data review groups to determine if analytical results are

acceptable. If recoveries of all quality control samples are within the method specified tolerances then

the data will be presented to laboratory project managers for final review. If recoveries of quality control

samples exceed the specified tolerances affected samples will be submitted for re-analysis.

Final review of analytical results will consist of the Laboratory Project Director's determination that all

analytical results of a sample(s) are consistent. If so, the data will be presented in a final report. If

discrepancies or deficiencies exist in the analytical results, corrective action will be taken. Audits of

final reports by the Laboratory Quality Assurance Officer may be conducted to determine the precision,

accuracy, completeness, and representativeness of sample analyses.

TOLEDO TIE TREATMENT SITE QUALITY ASSURANCE PROJECT PLAN FINAL PAGE 20 of 33

#### C.8.3 Data Reporting

Data reporting will be in accordance with the appropriate U.S.EPA method used for analysis. Laboratory reports shall, at a minimum, include the following:

- 1. narrative including statement of samples received, description and rationale for any deviations from approved methods/SOPs, summary of data quality, and documentation of any significant problems encountered during analysis;
- 2. documentation of laboratory events including dates of sample receipt, sample extraction, and sample analysis;
- 3. analytical data including results, detection limits, dilutions, etc.;
- 4. a summary of QA/QC results and supporting documentation;
- 5. a copy of the signed COC for samples submitted for analysis.

Laboratory reports should be signed by the laboratory's QAO and/or the laboratory director prior to being issued. Reports will be issued to HAI's QAO. Any draft reports should be clearly identified as such.

#### C.8.4 Data Validation

Analytical data will be reviewed according to the laboratory's data validation procedures outlined in Attachment A. After passing internal data validation, the data will be reported to HAI's QAO. Data will be reviewed by HAI's QAO to determine that proper preservation, holding times, and sample analysis procedures have been followed and are clearly documented. Additionally, the analytical results will be reviewed and compared to previous data, if any. Any questions regarding the data reports will be brought to the laboratory project manager's attention.

If additional data validation is required by the U.S.EPA or appropriate state agencies, the data will be submitted for validation to an independent third-party data validator or the U.S. EPA. The analytical laboratory is required to address any comments and correct any deficiencies identified in the data validation report.

HULL & ASSOCIATES, INC. TOLEDO, CHIO

TOLEDO TIE TREATMENT SITE QUALITY ASSURANCE PROJECT PLAN

FINAL PAGE 21 of 33

**C.9.0 INTERNAL QUALITY CONTROL** 

C.9.1 General

The purpose of internal quality control measures is to document the validity of analytical data generated

by the laboratory. Laboratory internal quality control may include, but is not limited to, the analysis of

method blanks, reference standards, analytical spikes, and surrogate spikes. Every analytical series will

include some of these controls depending on the analytical methods used. The internal quality controls

used by the laboratory will be combined so they are completely representative of the analytical task from

sample preparation and sample analysis.

The following sections present a summary of, and suggested frequencies for, various quality control

measures that may be used dependent upon the analytical method(s) selected. The laboratory's QAPP,

located in Attachment A, presents the actual quality control measures and frequencies that will be

employed by the laboratory.

C.9.2 Blank Samples

Blanks are used to assess contamination introduced in transit, storage, or in the laboratory. The types

and frequencies of laboratory blank samples are specified by the U.S.EPA methods used for analysis.

C.9.2.1 Method Blanks

Method blanks identify sources of contamination throughout the analytical process, whether a contribution

of specific analytes or a source of interference, which will need to be identified, isolated, and corrected.

To accomplish this, the method blank must be initiated at the beginning of the analytical process and

include all aspects of the analytical work. This includes all glassware, reagents, and instrumentation,

as well as any other possible source of contamination. Minimum method blank analyses will be one

method blank per analytical series at a frequency of one per 20 samples.

HULL & ASSOCIATES, INC. TOLEDO, CHIO

TOLEDO TIE TREATMENT SITE QUALITY ASSURANCE PROJECT PLAN FINAL. PAGE 22 of 33

C.9.2.2 Container Blank

The same concept for the method blank will apply to the sample bottles furnished from the laboratory. This container blank will be analyzed for each type of sample container as it would be used for collection. The frequency of analysis will extend to each lot of processed sample containers. At a minimum, the analysis of a container blank should be performed whenever the preparation process, preservation reagent, or the type of container changes.

C.9.2.3 Holding Blanks

Another type of method blank is a holding blank. Holding blanks are associated with volatile organic analyses and indicate possible cross-contamination among samples while stored at the laboratory. At least one holding blank, per each group of samples, will be generated and analyzed with the samples.

C.9.3 Reference Standards

œ۱

Reference standards are standards of known concentration and independent in origin from the calibration standards. These reference standards are generally available through the U.S.EPA, the National Bureau of Standards, or are specified by analytical methodologies. The purpose of a reference standard is to assess analytical proficiency within an analytical series including the preparation of calibration standards, the validity of calibration, sample preparation, instrument set-up, and the premises inherent in quantitation. Reference standards will be used in every analytical series except GC/MS and specific GC analyses, for which there are no reference standards.

A control chart will be maintained for analytes in which reference standards are used in their analyses. When a reference standard value exceeds the established warning limits, careful scrutiny will be given to the operations system, standards preparation, and procedures that were used in obtaining the result. If the value of the reference standard exceeds the established control limits, then the sample analysis will be stopped and corrective action will be initiated. Samples analyzed since the last passing reference standard will be re-analyzed following instrument recalibration.

HULL & ASSOCIATES, INC. APRIL 1998 PWM001D 001

TOLEDO, OHIO

TOLEDO TIE TREATMENT SITE
QUALITY ASSURANCE PROJECT PLAN
FINAL

PAGE 23 of 33

**C.9.4** Analytical Spikes

The purpose of an analytical spike is to assess the efficiency and proficiency of an analytical series. This

includes quantitation standards, sample preparation, instrument set-up, and the premises inherent in

quantitation. This control reflects the competency of sample analysis within an analytical series while

it is less sensitive in reflecting the conditions which are within the control of the analyst. The types and

frequencies of analytical spikes are specified by the U.S.EPA methods used for analysis.

C.9.4.1 Matrix Spike

Within an analytical series, a representative sample portion is designated as a separate sample and spiked

with known concentrations of the analytes under consideration. Advantages of spikes are that the spiked

portion is handled and prepared in exactly the same way as the samples. Sample related interference

affecting analysis will be reflected in the results from the spiked sample. Results of spikes exceeding

tolerances specified by the methods need to be evaluated thoroughly in conjunction with other measures

of control.

C.9.4.2 Surrogate Spike

Surrogates, which have properties similar to the analytes of interest, are compounds unlikely to be found

in nature. The intent of a surrogate spike is to provide broader insight to the proficiency and efficiency

of an analytical method on a sample specific basis. This control reflects analytical conditions which may

not be attributable to the sample matrix. If results of a surrogate spike analysis exceed method specified

tolerances, then the analytical results need to be evaluated thoroughly in conjunction with other control

measures. Re-analysis of the sample with additional controls, or different analytical methodologies, will

be necessary.

C.9.5 Replicate Analysis

Replicate analysis is a measure of analytical precision and can be limited in its scope. If used in

conjunction with reference standards or analytical spikes, it can measure the reliability of the analytical

systems. Replicate analyses can be significant in the interpretation of analytical results for samples with

complex matrices.

HULL & ASSOCIATES, INC. TOLEDO, OHIO

TOLEDO TIE TREATMENT SITE QUALITY ASSURANCE PROJECT PLAN FINAL PAGE 24 of 33

#### C.9.6 Calibration Check Standards

The purpose of a calibration check standard is to assess an instrument's stability. A calibration check standard will be analyzed at the beginning and end of an analytical series or periodically throughout large series of samples. Calibration check standard will be run after every ten samples. In analyses where internal standards are used, a calibration check standard need only be run at the beginning of an analytical series. If results of the calibration check standard exceed method specified tolerances, then samples analyzed since the last acceptable calibration check standard will be re-analyzed.

#### C.9.7 Internal Standards

Internal standards will be monitored when required by the method (e.g., U.S.EPA Method 624). The internal standard is present in all acquisitions with the exception of performance standards. The response of each compound within the internal standard is plotted on a control chart. The area of any compound cannot fall below 50% of its value in the preceding check standard, nor can it rise above 100% of its value. If internal standard areas in one or more samples exceed the specified tolerances, then the instrument will be recalibrated and all affected samples re-analyzed.

FINAL PAGE 25 of 33

C.10.0 PERFORMANCE AND SYSTEM AUDITS

C.10.1 General

Performance and system audits of both field and laboratory activities will be conducted to verify that

sampling and analysis are performed within the constraints of this plan. These audits can either be

conducted internally by field or laboratory staff or externally by state or federal agencies. The laboratory

will participate in any performance or system audit conducted or requested by HAI, appropriate state

agencies, or the U.S.EPA.

**C.10.2 Performance Audits** 

Performance audits will be conducted periodically to determine the accuracy of the total measurement

system(s) or components. In this program, blind performance evaluation samples, submitted by state and

federal agencies, are analyzed and evaluated throughout the year as part of an ongoing participation in

their certification programs. Any deficiencies in the results of these analyses are reported to the

laboratory and corrective action is initiated.

In addition to blind sample analyses, the laboratory will also participate in any audits from state and

federal agencies. These agencies submit a report noting any deficiencies and necessary corrective action.

The laboratory will respond with evidence of compliance within a limited time.

The laboratory also maintains a schedule of internal audits whereby each section of the laboratory is

audited by the Laboratory Quality Assurance Officer. When the audit is completed, a formal report will

be issued to the Laboratory Director. This report shall note any deficiencies and a follow-up date to

confirm corrective action.

C.10.3 System Audits

A system audit is an evaluation of the various components of the measurement system to asses their

proper selection and use. This includes a careful evaluation of all laboratory quality control measures.

System audits will be conducted internally by the laboratory.

HULL & ASSOCIATES, INC. TOLEDO, OHIO

TOLEDO TIE TREATMENT SITE QUALITY ASSURANCE PROJECT PLAN FINAL PAGE 26 of 33

#### C.10.4 Field Audits

Internal audits of field activities (sampling and measurements) will be conducted by HAI's QAO and/or FOC. These audits will include examination of field sampling records, field instrument operating records, sample collection, shipping and handling, COC, etc. These audits will occur at the onset of the project to verify that the established procedures are followed. Follow-up audits will be conducted to correct deficiencies, and to verify the QA/QC procedures are being maintained throughout the project. When an audit is completed a written report will be submitted to the PM.

HAI personnel will participate in any external audit requested by state and federal agencies. The results and recommendations or any external audit should be reported to HAI's QAO and/or PM in a timely manner so that corrective actions may be initiated.

TOLEDO TIE TREATMENT SITE QUALITY ASSURANCE PROJECT PLAN FINAL PAGE 27 of 33

#### **C.11.0 PREVENTIVE MAINTENANCE**

Field instruments and equipment will be maintained and serviced according to the manufacturer's instruction manual. A maintenance record for each instrument will be maintained. These records will include dates and descriptions of service and preventive maintenance. Major maintenance on any environmental instruments will only be performed by the manufacturer or trained technicians. Field personnel will be responsible for daily maintenance of all equipment in their possession. Critical spare parts (e.g., batteries/chargers, probes, screws, etc.) and tools will be kept with the equipment. Equipment problems will be reported to the QAO or the PM immediately.

Laboratory instruments will be maintained and serviced according to the individual instrument manuals. The laboratory's preventative maintenance procedures are documented in Attachment A.

HULL & ASSOCIATES, INC. TOLEDO, OHIO

#### C.12.0 DATA PRECISION, ACCURACY, AND COMPLETENESS

#### C.12.1 Precision

Precision is a measure of agreement between repetitive measurements under identical conditions. The overall precision of measurement data is a mixture of sampling and analytical factors. Analytical precision is much easier to control and quantify than sampling precision. Sampling precision may be determined by analyzing duplicates or replicate field samples and then creating and analyzing laboratory replicates from one or more of the field samples. The analytical results from the field duplicates or replicates provide data on overall measurement precision. Analytical results from the laboratory replicates provide data on analytical precision. Sampling precision can be calculated by subtracting the analytical precision from the overall measurement precision. For organic analyses, precision is reported as the Relative Percent Difference (%RPD) between matrix spike and matrix spike duplicate analysis. For metal analyses, precision is reported as %RPD between two duplicate samples. Acceptable limits for precision are specified by the U.S.EPA method. In the presence of outliers, corrective action will be taken including repairing instruments and/or re-analysis of the affected sample or samples. The following equation will be used in calculating %RPD:

$$%RPD = \underline{D_1 - D_2}_{X 100}$$

$$(D_1 + D_2)/2$$

RPD = Relative Percent Difference

D1 = First Sample Value

D2 = Second Sample Value (duplicate)

#### C.12.2 Accuracy

Accuracy is the difference between a measured value and the actual value, or the bias in a measurement system. Accuracy is difficult to measure for the entire data collection activity. Sources of error are the sampling process, field contamination, preservation, handling, sample matrix, sample preparation, and analysis techniques. Sampling accuracy can be assessed by evaluating the results of field and trip blanks.

TOLEDO TIE TREATMENT SITE QUALITY ASSURANCE PROJECT PLAN FINAL

PAGE 29 of 33

For the analytical procedures to be used, accuracy is assessed by Matrix Spike Percent Recovery. Acceptable limits for accuracy are specified by the U.S.EPA method. In the presence of outliers,

corrective action will be taken including fixing instruments and/or re-analysis of the affected sample or

samples. Matrix spike percent recovery will be calculated by the following equation:

 $MSPR = \underline{SSR-SR} \times 100$  SA

MSPR = Matrix Spike Percent Recovery

SSR = Spike Sample Results

SR = Sample Results

SA = Spike Added (concentration)

C.12.3 Completeness

Completeness will be reported as the percentage of the measurements judged to be valid. Completeness goals for this project are presented in Section C.3.5. Completeness will be calculated by the following

equation:

Completeness = 100 x (valid measurements / total measurements)

HULL & ASSOCIATES, INC TOLEDO, O'HO

PAGE 30 of 33

**C.13.0 CORRECTIVE ACTION** 

C.13.1 General

Corrective actions may be required for either analytical and equipment problems or noncompliance

problems. Analytical and equipment problems may occur during sampling and sample handling, sample

preparation, laboratory analysis, and data review. Noncompliance problems are associated with

nonconformance to this plan or the U.S.EPA methods being used.

C.13.2 Laboratory Corrective Action

When deficiencies or "out-of-control" situations exist, the laboratory will provide a means of detecting

and correcting these situations. An "out-of-control" situation is defined as data exceeding control limits.

Samples analyzed during "out-of-control" situations will be re-analyzed prior to reporting results. The

laboratory's corrective action procedures are documented in Attachment A. In general, there are several

levels of "out-of-control" situations that may occur in the laboratory during analysis.

C.13.2.1 Bench Level

Corrective action procedures will often be handled at the bench level. If an analyst finds a non-linear

response during calibration of an instrument, then the instrument will be recalibrated before sample

analysis. The problem may be corrected by a careful examination of the preparation or extraction

procedure, spike and calibration mixes, or instrument sensitivity. If the problem persists, it will be

brought to the management level.

C.13.2.2 Management Level

If resolution at the bench level was not achieved, or a deficiency is detected after the data has left the

bench level, then corrective action becomes the responsibility of the Laboratory Manager or Director.

Unacceptable matrix or surrogate spike recoveries detected by data review will be reported to the

Laboratory Manager. A decision to re-analyze the sample or report results will be made depending on

the circumstance.

HULL & ASSOCIATES, INC.

APRIL 1998
TOLEDO, OHIO

PWM001D.001

TOLEDO TIE TREATMENT SITE QUALITY ASSURANCE PROJECT PLAN FINAL PAGE 31 of 33

#### C.13.2.3 Receiving Level

If discrepancies exist in either the documentation of a sample or its container, a decision will be made after consulting with the appropriate management personnel. Decisions will be fully documented. Some examples of container discrepancies are broken samples, inappropriate containers, or improper preservation. In these cases, corrective action will involve the Laboratory Project Manager contacting HAI's QAO.

#### **C.13.3 Field Corrective Action**

Corrective actions for field equipment problems will consist of reporting the problem to the PM and/or the QAO so that maintenance can be performed or new equipment can be acquired. Noncompliance problems will be reported immediately to the QAO. The QAO will consult with the PC and corrective actions will be initiated. When warranted, the PM will report the nature of the noncompliance and corrective actions implemented to the appropriate state agencies and/or the U.S. EPA. The nature, extent, and corrective action for all noncompliances will be documented.

#### C.14.0 QUALITY ASSURANCE REPORTS TO MANAGEMENT

#### C.14.1 Internal Reporting

The Laboratory Quality Assurance Officer will report the status of the laboratory QA/QC program to the laboratory management. Each report will include:

- 1. periodic assessment of measurement data accuracy, precision, and completeness;
- 2. results of audits;
- 3. significant QA/QC problems and recommended solutions; and
- 4. resolutions of previously stated problems.

The laboratory will determine the content and frequency of these reports in accordance with its QAPP, which is included as Attachment A, and its SOPs. The laboratory will report to HAI's QAO when the results of HAI's samples have been affected by internal quality issues.

#### C.14.2 Additional Reporting

Laboratory analytical reports will include a summary of the quality assurance activities and quality control data for the project as related to sample analysis. The Laboratory Project Manager will report suspected field QA/QC problems to HAI's QAO.

HAI's QAO will report to HAI's PM when appropriate. These reports may be either oral or written depending upon the nature and complexity of the issues in the report.

#### C.15.0 REFERENCES

A variety of technical manuals, administrative documents, and publications were referred to in preparing this document. Some of the references consulted are presented below. Referenced documents and publications may or may not have been reviewed in their entirety. The guidelines and procedures presented in the documents and publications referenced have not been strictly adhered to unless stated otherwise.

- U.S.EPA. <u>Interim Guidelines and Specifications for Preparing Quality Assurance Project Plans</u>. EPA/600/4-83-004. February 1983.
- U.S.EPA. <u>Data Quality Objectives for Remedial Response Activities: Development Process.</u> EPA/540/6-87/003. March 1987.
- U.S.EPA. <u>Data Quality Objectives for Remedial Response Activities: Example Scenario</u>. EPA/540/6-87/004. March 1987.
- U.S.EPA. <u>A Compendium of Superfund Field Operations Methods</u>. EPA/540/P-87/001. December 1987.
- U.S.EPA. <u>Test Methods for Evaluating Solid Waste, Physical/Chemical Methods</u>. SW-846, 3rd Edition. September 1986.
- U.S.EPA. Methods for Chemical Analysis of Water and Wastes. EPA/600/4-79-020. March 1983.
- U.S.EPA. Quality Assurance/ Quality Control Guidance for Removal Activities. EPA/540/G-90/004. April 1990.

# Table C.1

U.S. EPA Method	<b>Parameters</b>
8260 (VOCs)	1,1-Dichloroethene
	Chloromethane
	Bromomethane
	Vinyl Chloride
	Chloroethane
	Acrolein
	Acrylonitrile
	Methylene Chloride
	Trichlorofluoromethane
	1,1-Dichloroethane
	Chloroform
	1,2-Dichloroethane
	1,1,1-Trichloroethane
	Dichlorodifluoromethane
	Methyl Iodide
	Acetonitrile
	Acetone
	Carbon Disulfide
	Propionitrile
	Allyl Chloride
	Methacrylonitrile
	2-Butanone
	Dibromomethane
	1,4-Dioxane
	trans-1,2-Dichloroethene
	cis-1,2-Dichloroethene
	Trichloroethene
	Benzene
	Toluene
	Chlorobenzene
	Xylene (total)
	Carbon Tetrachloride
	Bromodichloromethane
	1,1,2,2-Tetrachloroethane

# Table C.1

U.S. EPA Method	Parameters Parameters
8260 (VOCs) cont.	1,2-Dichloropropane trans-1,3-Dichloropropene Dibromochloromethane 1,1,2-Trichloroethane cis-1,3-Dichloropropene Bromoform Tetrachloroethene Ethylbenzene Isobutyl Alcohol Vinyl Acetate 2-Chloro-1,3-butadiene 1,2-Dibromoethane Methyl Methacrylate 1,1,1,2-Tetrachloroethane tans-1,4-Dichloro-2butene 1,2,3-Trichloropropane 2-Hexanone 4-Methyl-2-pentanone Ethyl Methacrylate Pentachloroethane 1,2-Dibromo-3 chloropropane Styrene
8270 (SVOCs)	aniline acetophenone bis(2chloroisopropyl)ether ethyl methanesulfonate methyl methanesulfonate N-nitrosodiethylamine N-nitrosomethylethylamine N-nitrosomorpholine N-nitrosopiperidine N-nitrosopyrolidine 2-picoline

# Table C.1

U.S. EPA Method	<u>Parameters</u>
8270 (SVOCs) cont.	o-toludine
	benzyl alcohol
	bis (2-chloroethoxy) methane
	bis (2-chloroethyl) ether
	2-chlorophenol
	2-methylphenol
	3- or 4-methylphenol
	1,2-dichlorobenzene
	1,3-dichlorobenzene
	1,4-dichlorobenzene
	hexachloroethane
	isophorone
	nitrobenzene
	2-nitrophenol
	N-nitrosodimethylamine
	N-nitrosodi-n-propylamine
	phenol
	pyridine
	2,6-dichlorophenol
	1,3-dinitrobenzene
	hexachloropropene
	isosafrole
	1,4-naphthoquinone
	N-nitrosodi-n-butylamine
	1,4-phenylenediamine
	safrole
	1,2,4,5-tetrachlorobenzene
	0,0,0-
	triethylphosphorothioate
	a,a-dimethylphenethylamine
	4-chloroaniline
	4-chloro-3-methylphenol
	2-chloronaphthalene
	2,4-dichlorophenol
	2,4-dimethylphenol
	dimethyl phthalate

#### Table C.1

U.S. EPA Method	<u>Parameters</u>
8270 (SVOCs) cont.	2,6-dinitrotoluene
	hexachlorobutadiene
	hexachlorocyclopentadiene
	2-methylnaphthalene
	naphthalene
	2-nitroaniline
	1,2,4-trichlorobenzene
	2,4,5-trichlorophenol
	2,4,6-trichlorophenol
	2,3,4,6-Tetrachlorophenol
	dimethoate
	1-naphthylamine
	2-naphthylamine
	5-nitro-o-toluidine
	pentachlorobenzene
	phenacetin
	tetraethyldithiopyrophosphate
	1,3,5-trinitrobenzene
	diallate (trans/cis)
	acenaphthene
	acenaphthylene
	4-bromophenyl phenyl ether)
	4-chlorophenyl phenyl ether
	dibenzofuran
	diethyl phthalate
	4,6-dinitro-2-methylphenol
	2,4-dinitrophenol
	2,4-dinitrotoluene
	fluorene
	hexachlorobenzene
	3-nitroaniline
	4-nitroaniline
	4-nitrophenol
	N-nitrosodiphenylamine
	2-acetylaminofluorene

# Table C.1

# Parameter List

U.S. EPA Method	Parameters
8270 (SVOCs) cont.	
8270 (SVOCs) cont.	4-aminobiphenyl chlorobenzilate p-(dimethylamino)azobenzene 7,12-dimethylbenz(a)anthracene 3,3'-dimethylbenzidine isodrin 3-methylcholanthrene 4-nitroquinoline-1-oxide pentachloronitrobenzene pronamide methapyrilene thionazin anthracene benzo (a) anthracene benzo (b) fluoranthene benzo (k) fluoranthene benzo (ghi) perylene benzo (a) pyrene bis (2-ethylhexyl) phthalate chrysene di-n-butyl phthalate dibenz (a,h) anthracene 3,3'-dichlorobenzidine di-n-octyl phthalate fluoranthene indeno (1,2,3-cd) pyrene
8081(Organochlorine Pesticides)	pentachlorophenol phenanthrene pyrene  Kepone Heptachlor Aldrin Dieldrin

HULL & ASSOCIATES, INC. TOLEDO, OHIO

APRIL 1998 PWM001D.001

# Table C.1

.S. EPA Method	Parameters
8081 (Organochlorine Pesticides) cont.	Endrin
( - 5	DDT
	Alpha BHC
	Chlordane
	Beta BHC
	Gamma BHC - Lindane
	Delta BHC
	Heptachlor Epoxide
	DDE
	DDD
	Methoxychlor
	Toxaphene
	Endosulfan I
	Endosulfan II
	Endosulfan Sulfate
	Endrin Aldehyde
	PCB-1016
	PCB-1221
	PCB-1232
	PCB-1242
	PCB-1248
	PCB-1254
	PCB-1260
8141 (Organochlorine Pesticides)	Ethyl Parathion
	Famphur
	Phorate
	Disulfoton
	Methyl Parathion
8252 (Chlorinated Herbicides)	2,4-D
The committee of the control of the	2,4,5-TP
	2,4,5-T
	Dinoseb
	Hexachlorophene
	Hexacinorophene

#### Table C.1

Parameter List				
U.S. EPA Method	Parameters			
7470 (Inorganics/Metals)	Arsenic			
( <b>g</b> )	Barium			
	Cadmium			
	Chromium			
	Copper			
	Lead			
	Mercury			
	Selenium			
	Silver			
	Zinc			

Table C.2

Sample Containers and Preservation

Parameter	Containers	Preservatives/ Max Holding Time
SVOCs (8270)	Soils - one 8oz. wide mouth jar.	Cool 4°C / 10 (40) days <sup>1</sup>
	Aqueous - three 1 liter bottles.	Cool 4°C / 5 (40) days <sup>1</sup>
Pesticides (8081)	Soils - one 8oz. wide mouth jar.	Cool 4°C / 10 (40) days <sup>1</sup>
	Aqueous - three 1 liter bottles.	Cool 4°C / 5 (40) days <sup>1</sup>
Metals (6000/7000)	Soils - one 8 oz. wide mouth jar	Cool 4°C
	Aqueous - 250 ml plastic	$HNO_3 pH > 2$
TCLP (1311)	Soils - two 8 oz. wide mouth jar	Cool 4°C
· ·	Aqueous - three 1 liter bottles.	Cool 4°C
Reactive sulfide	Soils - one 4 oz. wide mouth jar	Cool 4°C
	Aqueous - 250 ml plastic	Cool 4°C
Reactive cyanide	Soils - one 4 oz. wide mouth jar	Cool 4°C
·	Aqueous - 250 ml plastic	Cool 4°C
pH (150.1)	Soils - one 2 oz. wide mouth jar	Cool 4°C
	Aqueous - 250 ml plastic	Cool 4°C
Flash point (1010)	Soils - one 2 oz. wide mouth jar	Cool 4°C
	Aqueous - 250 ml plastic	Cool 4°C

Note:

All holding times are from time of sample collection. This list represents typical sample containers that may be supplied for this project. The contracted laboratory will provide a detailed sheet describing the types and number of containers sent for each analysis with each sample kit. The actual number of containers may be less as several parameters may be combined into a single container. The laboratory will also provide any required preservatives and instructions for preservation. If pre-preserved bottles are supplied, they will be clearly identified on the sampling container.

<sup>&</sup>lt;sup>1</sup>5 or 10 days pre-extraction / 40 days post extraction.

#### Table C.3

#### Sample Types and Codes

Code **Type** Water Monitoring Well MW(designated monitoring well 1.D.) Surface Water SW(designated monitoring station)) Sediment **SED** Soil Surface Soil SS(designated sampling location 1.D.) Soil Boring SB(designated boring I.D.) Drum Sample DS QA/QC **Duplicate** Same as well code FB(number)1 Field Blank Trip Blank TB(number)1

Numbers for field and trip blanks begin at 1 and increase by increments of 1 every time a new blank is collected. The number will begin at 1 at the beginning of every day.

#### Table C.4

#### **Field Equipment List**


- 1 Visqueen
- 2. non-phosphatic soap
- 3. Distilled water
- 4. Trash bags
- 5. Analytical containers
- 6. Shipping containers/coolers
- 7. Ice
- 8. Indelible ink pens
- 9. Clear packing tape
- 10. Health and safety equipment
- 11. Decon supplies
- 12. Tools

Note: This represents a general list of sampling equipment required for this project. Additional equipment required for specific tasks is presented in the FSAP.

#### Table C.5

# **Analytical Procedures**

U.S. EPA Methods	Parameters
1311 - 8240	TCLP volatiles
1311 - 8270	TCLP semi-volatiles
1311 - 8080 & 8150	TCLP pesticides/herbicides
1311 - 6000/7000 series	TCLP metals
1010	flash point
150.1	рН
SW-846 Section 7.3.3.2	reactive cyanide
SW-846 Section 7.3.4.1	reactive sulfide
6000/7000 series	metals
8260	volatile organics
8270	semi volatile organics
8081	organochlorine pesticides/PCBs
8141	organophosphorous pesticides
8252	chlorinated herbicides
9012	cyanide
9030	sulfide



**ATTACHMENT À** 

LANCASTER LABORATORIES' QUALITY ASSURANCE PLAN

HULL & ASSOCIATES, INC. TOLEDO, OHIO

QUALITY ASSURANCE PROJECT PLAN PWM001D.001

#### LABORATORY QUALITY ASSURANCE PLAN

DECEMBER 12, 1991 REVISED: May 12, 1997

**WARNING:** The information contained herein is of a highly confidential and proprietary nature. Lancaster Laboratories specifically prohibits the dissemination or transfer of this information to any person or organization not directly affiliated with the project for which it was prepared.



Section No. 1 Revision No. 4 Date: 05/12/97 Page 1 of 1

# 1. Laboratory Quality Assurance Plan

This document provides the laboratory portion of the response to EPA's *Interim Guidelines and Specifications for Preparing Quality Assurance Project Plans* QAMS-005/80, Sections 5.1 through 5.16 as revised December 29, 1980, and EPA-600/4-83-004, February 1983. Guidance was also obtained from *Preparation Aids for the Development of Category 1 Quality Assurance Project Plans*, Office of Research and Development, USEPA, EPA/600/8-91/003, February 1991.

As much as possible, the procedures in this document have been standardized to make them applicable to all types of environmental monitoring and measurement projects. However, under certain site-specific conditions, all of the procedures discussed in this document may not be appropriate. In such cases it will be necessary to adapt the procedures to the specific conditions of the investigation.

Quality Assurance Officer. Johnson 7. Journ

Section No. 2 Revision No. 5 Date: 05/12/97 Page 1 of 1

		<u>Section</u>	<b>Pages</b>	Revision	<u>Date</u>
	1.	Title Page	1	4	05/12/97
	2.	Table of Contents	1	5	05/12/97
***	3.	Project Description	1	3	05/12/97
<b>41</b>	4.	Project Organization and Responsibility	4	4	05/12/97
**	5.	QA Objectives for Measurement Data, in terms of precision, accuracy, completeness, representativeness, and comparability	3	3	05/12/97
4	6.	Sampling Procedures	2	2	07/25/95
	7.	Sample Custody	26	4	05/12/97
484	8.	Calibration Procedures and Frequency	5	2	07/25/95
#	9.	Analytical Procedures	9	5	05/12/97
	10.	Data Reduction, Validation, and Reporting	6	2	07/25/95
	11.	Internal Quality Control Checks	9	3	05/12/97
•	12.	Performance and Systems Audits	18	3	10/23/96
	13.	Preventive Maintenance	3	2	07/25/95
•	14.	Specific Routine Procedures Used to Assess Data	3	2	07/25/95
•	15.	Precision, Accuracy, and Completeness Corrective Action	3	2	07/25/95
	16.	Quality Assurance Reports to Management	1	2	07/25/95
•	Appe	ndix A - CLP Forms Inorganics and Organics	70	2	07/25/95

Section No. 3 Revision No. 3 Date: 05/12/97 Page 1 of 1

#### 3. Project Description

Tests will be performed according to the analytical methodology set forth in the USEPA Contract Laboratory Program Statement of Work\*. The USEPA-CLP-SOW provides specific analytical procedures to be used and defines the specific application of these procedures. Proven instruments and techniques will be used to identify and measure the concentrations of volatiles, semivolatiles, and pesticide compounds and/or the inorganic elements. The laboratory will employ state-of-the-art GC/MS and/or GC procedures to perform all organic analyses, including all necessary preparation for analysis. Inorganic analyses will be performed using graphite furnace atomic absorption spectrophotometry (AA), inductively coupled plasma spectroscopy, cold vapor AA, or flame AA. The client is responsible for providing specifics on the project site.

\*USEPA-CLP-SOW for Organics, Document No. OLM03.2, USEPA-CLP-SOW for Inorganics, Document No. ILM04.0, or most recent revision unless otherwise requested by the client.

Section No. 4 Revision No. 4 Date: 05/12/97 Page 1 of 4

#### 4. Project Organization

The objectives of the laboratory Quality Assurance Program are to establish procedures which will ensure that data generated in the laboratory are within acceptable limits of accuracy and precision, to ensure that quality control measures are being carried out, and to ensure accountability of the data through sample and data management procedures. To this end, a Quality Assurance Department has been established. The Quality Assurance Officer reports directly to the President of Lancaster Laboratories and has no direct responsibilities for data production, thus avoiding any conflict of interest.

The attached organizational charts show key managerial personnel. Resumes of key individuals may be found in the enclosed *Qualifications Manual*.


The Sample Administration Group will be responsible for receiving samples, signing the external chain of custody, checking sample condition, assigning unique laboratory sample identification numbers, and initiating internal chain-of-custody forms. Sample Support personnel will be responsible for assigning storage locations, checking and adjusting preservation, homogenizing the sample as needed, and sample discard.

Group leaders listed in each technical area are responsible for performing laboratory analyses, quality control as specified in the methods, instrument calibration, and technical data review. Data is reported using a computerized sample management system, which tracks sample progress through the laboratory and generates client reports when all analyses are complete. Quality control data is entered onto the same system for purposes of charting and monitoring data quality.

The Quality Assurance Department is responsible for reviewing quality control data, conducting audits in the laboratory and reporting findings to management, maintaining current copies of all analytical methods, maintaining copies of computer code used to calculate and report results, submitting blind samples to the laboratory, and ensuring that appropriate corrective action is taken when quality problems are observed.


Section No. 4 Revision No. 4 Date: 05/12/97 Page 2 of 4

Data package deliverables are available upon request. The Quality Assurance Department reviews the contents of the deliverables for completeness and to be sure that all quality control checks were performed and met specifications. This step includes review of holding times, calibrations, instrument tuning, blank results, duplicate results, matrix spike results, and surrogate results. Every attempt to meet specifications will be made, and any item outside of the specifications will be noted in the narrative. The laboratory will not validate data with regard to usability since this generally requires specific knowledge about the site.



Section No. 4
Revision No. 4
Date: 05/12/9
Page 3 of 4

#### **Environmental Sciences**



Section No. 4
Revision No. 4
Date: 05/12/97
Page 4 of 4

Section No. 5 Revision No. 3 Date: 05/12/97 Page 1 of 3

#### 5. QA Objectives for Measurement Data

Quality assurance is the overall program for assuring reliability of monitoring and measurement data. Quality control is the routine application of procedures for obtaining set standards of performance in the monitoring and measurement process. Data quality requirements are based on the intended use of the data, the measurement process, and the availability of resources. The quality of all data generated and processed during this investigation will be assessed for precision, accuracy, representativeness, comparability, and completeness. These specifications will be met through precision and accuracy criteria as specified in Section 11. Detection limits are presented in Section 9.

<u>Precision</u> - Precision is determined by measuring the agreement among individual measurements of the same property, under similar conditions. The laboratory objective is to equal or exceed the precision demonstrated for the applied analytical method on comparable samples. The degree of agreement is expressed as the relative percent difference (RPD%). Evaluation of the RPD% is based on the criteria set forth in the Contract Laboratory Program (CLP) for organic and inorganic analyses. External evaluation of precision is accomplished by analysis of standard reference material and interlaboratory performance data.

Accuracy - Accuracy is a measure of the closeness of an individual measurement to the true or expected value. Analyzing a reference material of known concentration or reanalyzing a sample which has been spiked with a known concentration/amount is a way to determine accuracy. Accuracy is expressed as a percent recovery (%R). Evaluation of the %R is based on the criteria established for the CLP for organic and inorganic analyses.

Representativeness - Representativeness expresses the degree to which data accurately represents the media and conditions being measured. The representativeness of the data from the sampling site will depend on the sampling procedure. Sample collection is the responsibility of the client. Samples will be homogenized, if required, as part of the laboratory sample preparation. By comparing the quality control data for the samples against other data for similar samples analyzed at the same time, representativeness can be determined for this objective.

Section No. 5 Revision No. 3 Date: 05/12/97 Page 2 of 3

<u>Comparability</u> - Comparability conveys the confidence with which one set of data can be compared to another. The analytical results can be compared to other laboratories by using traceable standards and standard methodology and consistent reporting units. The Laboratory Quality Assurance Program documents internal performance, and the interlaboratory studies document performance compared to other laboratories.

<u>Completeness</u> - Completeness is a measure of the quantity of valid data acquired from a measurement process compared to the amount that was expected to be acquired under the measurement conditions. The completeness of an analysis can be documented by including in the data deliverables sufficient information to allow the data user to assess the quality of the results. Additional information will be stored in the laboratories archives, both hard copy and magnetic tape. Quality Assurance standard operating procedures (SOPs) are in place to provide traceability of all reported results.

To ensure consistent attainment of the quality assurance objectives, SOPs are in place detailing the requirements for the correct performance of laboratory procedures. The laboratory SOPs fall under five general categories:

- 1. Corporate policy
- 2. Quality assurance
- 3. Sample administration
- 4. General laboratory procedures
- 5. Analytical (i.e., methods, standard preps., instrumentation)

All SOPs are approved by the QA Department prior to implementation. The distribution of current SOPs and archiving of outdated ones are controlled through a master file. Table 5-1 provides an index of QA SOPs in place in support of the Quality Assurance objectives. These requirements are supplemented by the procedures in the laboratory and analytical SOPs.

Section No. 5 Revision No. 3 Date: 05/12/97

Page 3 of 3

Table 5-1					
Document #	Document # Document Title				
QA-101	Sample Collection				
QA-102	Sample Log-in				
QA-103	Sample Storage and Disposal				
QA-104	Internal Chain-of-Custody Documentation				
QA-105	Analytical Methods Manual				
QA-106	Validation and Authorization of Analytical Methods				
QA-107	Analytical Methods for Nonstandard Analyses				
QA-108	Subcontracting to Other Laboratories				
QA-109	Laboratory Notebooks, Logbooks, and Documentation				
QA-110	Reagents				
QA-111	Instrument and Equipment Calibration				
QA-112	Instrument and Equipment Maintenance				
QA-113	Data Entry and Verification				
QA-114	Data Storage and Security				
QA-115	Quality Control Records				
QA-116	Investigation and Corrective Action of Unacceptable Quality Control Data				
QA-117	Personnel Training Records and Curriculum Vitaes				
QA-118	Quality Assurance Audits				
QA-119	Proficiency Samples				
QA-120	Documentation of Programming for the Sample Management System				
QA-121	Guidelines for the Development, Validation, Implementation, and Maintenance of Computer Systems Used with CLP, GLP, and GMP Data				
QA-122	Investigation and Corrective Action Reporting for Laboratory Problems				

Section No. 6 Revision No. 2 Date: 07/25/95 Page 1 of 2

#### 6. Sampling Procedures

In order for meaningful analytical data to be produced, the samples analyzed must be representative of the system from which they are drawn. It is the responsibility of the client to ensure that the samples are collected according to accepted or standard sampling methods.

The laboratory will provide the appropriate sample containers, required preservative, chain-of-custody forms, shipping containers, labels, and seals. The majority of sample containers are purchased precleaned I-Chem™ Series 200 or equivalent. Any reused bottles are cleaned in-house following laboratory standard operating procedures. Special containers with traceability documentation are available upon request. Because the laboratory does not stock this type of container, 1-month prior notice is required.

Each lot of preservative will be documented and checked for contaminants before use. The appropriate bottle will be preserved with the new preservative and filled with deionized water to represent a sample. A similar container (that does not contain preservative) will be filled with deionized water to be used as a blank check. Analysis results are documented for each preservative lot number.

Trip blanks will be prepared by the laboratory and accompany sample containers at the project required frequency. Analyte free water will also be provided for field blanks.

A list of containers, preservatives, and holding times follows in Table 6-1.

Section No. 6 Revision No. 2 Date: 07/25/95

Page 2 of 2

Table 6-1							
	Sample Containers, Preservatives, and Holding Times for Aqueous and Solid Samples						
Fraction	Vol. Req (mL) Wt. Req. (g)	Container P = Plastic G = Glass	Preservation <sup>a</sup>	Holding From D Reco Water	ate of		
Volatiles	3 × 40 mL 100 g	G	Cool, 4°Cb pH <2 w/HCl	10 Da	10 ys		
Pesticides	2 × 1000 mL 100 g	G	Cool, 4°C <sup>b</sup>	5 10  Days to  extraction <sup>e</sup>			
Acid/Base Neutrals	3 x 1000 mL 100 g	G	Cool, 4°C <sup>b</sup>	5 Day extrac			
Metals	1000 mL 100 g	P,G	HNO <sub>3</sub> to pH <2	6 Mor Hg 26	-		
Cyanide	1000 mL 100 g	P,G	Cool, 4°C NaOH to pH >12	12 12 Days			

<sup>&</sup>lt;sup>a</sup>pH Adjustment with acid/base is performed on water samples only.

NOTE: For volatiles analysis, the container should be filled completely, with no headspace. All sample containers, preservatives, and mailers will be supplied at no additional charge upon request, except for the special containers with traceability documentation. There is an additional charge for this type of container.

<sup>&</sup>lt;sup>b</sup>Sodium thiosulfate needed for chlorinated water samples

<sup>&</sup>lt;sup>c</sup>Assuming delivery of samples is within 2 days of sampling.

<sup>&</sup>lt;sup>d</sup>Samples will be analyzed as soon as possible after receipt. The times listed as the maximum times that samples will be held before analysis and still be considered valid.

Extracts of either water or soil/sediment samples must be analyzed within 40 days following extraction.

Section No. 7 Revision No. 4 Date: 05/12/97 Page 1 of 26

#### 7. Sample Custody

Samples are unpacked and inspected in the sample receipt area. At this time, the samples are examined for breakage and agreement with the associated client paperwork. The cooler temperatures will be checked upon receipt and recorded. As the samples are unpacked, the sample label information will be compared to the chain-of-custody record and any discrepancies or missing information will be documented. If necessary, the cooler will be closed and placed in cold storage until instructions and resolution of any discrepancies are received from the client.

A member of our Sample Administration Group will act as sample custodian for the project. To ensure accountability of our results, a unique identification number is assigned to each sample as soon as possible after receipt at the laboratory. When samples requiring preservation by either acid or base are received at the laboratory, the pH will be measured and documented, with the exception of samples designated for volatile analysis. Samples requiring refrigeration will be stored in our walk-in cooler which is maintained at 2° - 4°C. The use of our computer system in tracking samples (by the Lancaster Labs sample number assignment) will control custody of the sample from receipt until the time of its disposal. The security system on our laboratory building allows us to designate the entire facility as a secure area since all exterior doors are either locked or attended. Therefore, hand-to-hand chain of custody is not part of our routine procedure, but is available upon request. If requested, hand-to-hand chain of custody will be provided as per attached SOP-QA-104, "Chain-of-Custody Documentation." The laboratory chain of custody will begin with the preparation of bottles. The procedures for sample log-in, storage, and chain-of-custody documentation are detailed in the QA standard operating procedures included in Section No. 7 (QA-102, QA-103, and QA-104). Examples of sample labels and a custody seal are shown in Figure 7.1.

Section No. 7 Revision No. 4 Date: 05/12/97

Page 2 of 26



Where quality is a science.

SOP-QA-102.01 Supersedes Date: 02/15/96 Effective Date: MAR 1 4 1997

Page 1 of 5

#### **QUALITY ASSURANCE OPERATIONS MANUAL**

#### Sample Log-In

#### Purpose:

In order to provide accountability of our results, protect client confidentiality, and to prevent sample loss/mix-up, a continuous and unique Lancaster Laboratories (LL) identification number is assigned to each sample upon laboratory receipt.

#### Scope:

This SOP will cover the procedure used to log client samples into the computerized sample management system (SMS) after receipt. The Sample Administration Group is responsible for laboratory sample log-in. Sample Administration has procedures to define this sample entry process.

This procedure applies only to samples which are logged into and tracked by the SMS. There are only a few cases where samples may not be tracked using the SMS. These include samples which will be stored for a long period of time prior to analysis, (e.g., stability storage) or for special project samples that could be reported in a narrative research and development style report instead of our usual analytical reports. Written procedures for tracking samples not entered into the SMS are developed by the technical department responsible for the project or analysis of those samples.

#### **Personnel Training and Qualifications:**

Training in sample log-in is performed in accordance with Sample Administration training procedures.

Section No. 7 Revision No. 4 Date: 05/12/97 Page 3 of 26

SOP-QA-102.01

Supersedes Date: 02/15/96 Effective Date: MAR 1 4 1997

Page 2 of 5

#### Procedure:

- 1. All samples received for testing shall be delivered to the Sample Administration Department immediately upon arrival at the laboratory.
- 2. The Sample Administration Department will be responsible for unpacking and organizing the samples.
- 3. Client correspondence relating to the group of samples shall also be transferred to the Sample Administration Department. This may include purchase orders, quotations, letters, phone logs, and Incoming Sample Activity Records (ISARs).
- 4. Personnel of the Sample Administration Group shall log the samples into the SMS as soon as practical after receipt. Samples awaiting log-in are stored in temporary holding areas, at required temperature, to maintain the sample integrity. At the time of entry the computer will assign a unique identification number to each sample. Samples can be received at the laboratory 7 days a week, 24 hours a day, 365 days of the year. Samples should be logged in on the same day as they are received with the following exceptions:
  - a. Samples received on a holiday will not be logged-in until the next normal work day. Samples received from 6 p.m. on Saturday through 11 p.m. on Sunday will be logged-in Sunday evening by third shift Sample Administration personnel.
  - b. Samples submitted by clients which do not identify the type of testing to be performed or with unclear or incomplete paperwork documentation -Every effort will be made to contact the client on the same day of sample receipt. In this situation, the samples will be tracked in a hold database. The group of samples will be assigned a hold number. This database is maintained by the Sample Administration Group.

If same day sample log-in is not possible, all specified and appropriate storage requirements will be observed (e.g., refrigeration).

Section No. 7 Revision No. 4 Date: 05/12/97 Page 4 of 26

SOP-QA-102.01

Supersedes Date: 02/15/96 Effective Date: MAR 1 4 1997 Page 3 of 5

- 5. Upon assignment of a sample number, the SMS will generate a label which shall be attached to the sample container. Every effort will be made as to not obscure the client label. The information on the sample label will include the LL sample number, the client name, the storage location, the analyses requested, a bottle code indicating container and preservative type, if applicable, a unique bar code (used for samples stored in the Automated Sample Retrieval and Storage System [ASRS]), and any applicable notes to laboratory personnel.
- 6. Preservation, homogenization, and subsampling, if necessary, will be the responsibility of the Sample Support Group, or the testing laboratory. SOPs are in place within the group to define these procedures. A list of preservatives required for routine environmental analyses may be found in the Environmental Schedule of Services. A preservation, sulfate, and chlorine check shall be performed immediately after sample log-in for all applicable environmental samples.
- 7. After all above steps are performed, as required, samples shall be stored in an assigned storage location or taken to the laboratory for testing.
- 8. The next working morning, after sample log-in, a copy of an entry acknowledgment will print from the SMS. The acknowledgment is a hard copy record of the sample entry. It will summarize, the LL sample number, the sample(s) submitted in an entry group, the test(s) to be performed, the client requesting the work, the account to be billed for the work, and the unique sample identifications assigned by the client. This acknowledgment is mailed to the client to confirm sample receipt and entry.
- 9. Another copy of the sample acknowledgment will print and be designated as the laboratory copy. This acknowledgment, in addition to client paperwork, will audited by three levels of personnel after the entry process:
  - a. Sample Administration will audit to ensure that the entry corresponds to client supplied paperwork and/or quotations.

Section No. 7 Revision No. 4 Date: 05/12/97 Page 5 of 26

SOP-QA-102.01

Supersedes Date: 02/15/96 Effective Date: MAR 1 4 1997

Page 4 of 5

- b. Client Services will audit to ensure the entry is reflective of client documentation and that additional client/project requirements were communicated and taken into consideration. They will also verify that account and billing information is accurate.
- c. Technical centers will assure appropriate preparation and analysis set-up steps have been added to the entry. They will also verify that project and technical requirements have been taken into consideration from a technical point of view.

Each reviewer will initial the top of the SA file copy of the acknowledgment to document their review. Additional copies of this acknowledgment can be made for laboratory personnel.

10. The LL sample number assigned to each sample shall be used to identify the sample in all laboratory records, including laboratory notebooks, instrument printouts, and laboratory final reports. The sample number will also be used to identify all additional containers of the sample which may be created during sample preparation and analysis. This will include subsamples, extracts, and digests.

#### **Revision Log:**

Initiated Date: 03/87

Ver. # Effective Date Change
00 02/15/96 Previous issue

Section No. 7 Revision No. 4 Date: 05/12/97 Page 6 of 26

SOP-QA-102.01

Supersedes Date: 02/15/96 Effective Date: MAR 1 4 1997

Page 5 of 5

Ver. # Effective Date Change

01 MAR 1 4 1997 Major changes

Major changes are as follows:

- · Expanded upon the scope of the procedure
- Added section about printing and auditing of the sample entry acknowledgment
- Added Personnel Training and Qualification section
- Removed specifics on how to document preservation checks

SOPQA102.DOC 031197

Prepared by: Kathy D. Wetsel Date: 3/11/97

Approved by: Karlern M. Loewen Date: 3/1/97

Approved by: 3, hila Hershen Date: 3/12/97

Section No. 7 Revision No. 4 Date: 05/12/97 Page 7 of 26



Where quality is a science.

Initiated Date: 03/87

Effective Date: 007 0 1 1996

# QUALITY ASSURANCE OPERATIONS MANUAL SOP-QA-103

Title: Sample Storage and Disposal

#### Purpose:

Sample integrity can be compromised by improper storage conditions. The objective of these procedures is to prevent samples from deteriorating prior to analysis. The computerized sample management system (CSMS) is used to assign storage locations and to monitor the orderly storage of samples in locations from which they are easily retrieved for analysis or discard at the appropriate date.

#### Scope:

This SOP will outline procedures used in storing samples, retrieving and returning samples for analysis, and discarding samples when their holding time expires.

#### Procedures:

1. Personnel from Sample Administration will designate the approximate size and type (e.g., refrigerator, freezer or room temperature) of sample storage required for each group of samples as they are logged onto the CSMS. The computer will assign the storage location and record the length of time the sample must be retained after the analysis report has been issued. Samples will be stored in the assigned location. If the location is not suitable (e.g., insufficient space), the storage location may be changed using the manual override on the computer. If refrigerated space has been requested and all the computerized refrigerator locations are occupied, samples will be assigned locations in overflow refrigerators and will be tracked using a manual system until computerized locations are available.

Section No. 7 Revision No. 4 Date: 05/12/97 Page 8 of 26

SOP-QA-103 Initiated Date: 03/87

Effective Date: 0CT 0 1 1996

Page 2 of 4

2. Analysts requiring the use of a sample may determine its location by referring to the daily sample status sheet. There are varying degrees of security on sample storage locations. The procedures for removal of samples from these locations are as follows:

- a. Free access locations are those which are neither locked nor attended by a sample custodian. These areas are usually located within an individual group's laboratory and samples may be removed from and returned to these locations without documentation. However, if the sample must be taken out of the laboratory, documentation may be requested. Care shall be exercised in returning the sample to its appropriate location.
- b. Controlled access areas are attended by a sample custodian and are usually large areas used by more than one group. Samples stored in controlled access areas can be removed only after requisitioning the sample via the CSMS. The sample custodian will retrieve the requisitioned samples from the storage locations and scan the bar code label. This process documents the sample transfer from the sample custodian to laboratory personnel. After use, the samples are returned to the sample storage center, scanned by the sample custodian and returned to the designated storage location. Only Sample Administration personnel shall be admitted to controlled access areas. The only exception to this rule will be during weekend hours when no sample custodians are on duty. During these hours, samples must be requisitioned as above, but analysts must retrieve the samples themselves by obtaining a key to the controlled access area from the security desk. Samples must be scanned out as above. After use, samples must be scanned in and placed on the return cart inside WK. Sample custodians will return these samples to their location when they come on duty.

Section No. 7 Revision No. 4 Date: 05/12/97 Page 9 of 26

SOP-QA-103 Initiated Date: 03/87

Effective Date: 0CT 0 1 1996

Page 3 of 4

- c. Locked storage areas are available in several individual lab areas. Access to these storage areas is limited to analysts who are responsible for the analysis of the samples stored there. These areas are locked when the laboratories are unattended; keys are available from members of the department where they are located. Samples are removed and returned as needed by analysts.
- d. Forensic storage areas are locked and admission to these areas is permitted only to sample custodians. Most of the samples stored in these areas require strict chain-of-custody documentation as outlined in SOP-QA-104, "Internal Chain-of-Custody Documentation," and should be requisitioned as described in b. above. Samples may not be removed or returned to these areas without signing chain-of-custody forms.
- To prevent unnecessary deterioration of the samples, the aliquots needed for analysis shall be removed and the sample returned to storage with a minimum of delay.
- 4. Sample Administration will generate a discard list of samples with retention dates that have expired. The retention dates are based upon client requirements or defaulted to a given number of days past the date when the report is generated, if no client requirements were given. These samples will be removed from storage by a member of Sample Support or a member of the department responsible for the given storage location. Hazardous samples shall either be returned to clients, decontaminated or disposed of at the direction of supervisory personnel. Other samples will be discarded or returned to the client, if requested. Prior to discarding each sample, the bar code will be scanned to prevent discard of the wrong sample.
- The temperature of each refrigerator or freezer used for storing samples or reagents requiring temperature control should be checked during each normal working day by an assigned member of the group responsible for the

Section No. 7 Revision No. 4 Date: 05/12/97 Page 10 of 26

SOP-QA-103

Initiated Date: 03/87

Effective Date: OCT 0 1 1996

Page 4 of 4

samples stored within and recorded on a log posted on the outside of the unit. Units containing samples requiring more complete documentation of storage conditions are monitored by use of a computerized recording device or a temperature wheel. Refrigerator temperatures should be maintained at 2° to 4°C and freezer temperatures should be maintained at -15° ± 5°C, unless otherwise specified in a client-supplied method or protocol. If the temperature recorded does not fall within these ranges, the Maintenance Department should be contacted. Any repairs should be recorded and filed with the temperature log. All documentation of temperature checks and maintenance shall be kept in ink and any changes made shall follow the error correction procedure given in SOP-QA-109, "Laboratory Notebooks and Documentation."

SOPQA103.DOC 091196

Prepared/by:

OSEWEN)

Date:

Section No. 7 Revision No. 4 Date: 05/12/97 Page 11 of 26



Where quality is a science.

SOP-QA-104.01

Supersedes Date: 12/01/95 Effective Date: MAR 1 4 1997 Page 1 of 16

. ago . o. . o

# QUALITY ASSURANCE OPERATIONS MANUAL Internal Chain-of-Custody Documentation

## Purpose:

In order to demonstrate reliability of data which may be used as evidence in a legal case or required by a regulatory agency or client, an accurate written record tracing the possession of samples must be maintained from the time they are received at the aboratory until the last requested analysis is verified. The chain of custody is to ensure traceability of samples while they are in the possession of the laboratory.

# Scope:

Procedures for initiating and maintaining chain-of-custody (COC) documentation are described in this procedure.

#### Definition:

A sample is in custody if it is in any one of the following states:

- 1. In actual physical possession.
- 2. In view after being in physical possession.
- 3. Locked up so no one can tamper with it.
- 4. In a secured area, restricted to authorized personnel (e.g., in the ASRS system).

Section No. 7 Revision No. 4 Date: 05/12/97 Page 12 of 26

SOP-QA-104.01

Supersedes Date: 12/01/95 Effective Date: MAR 1 4 1997

Page 2 of 16

#### A. Procedure

- 1. Chain-of-custody documentation shall be kept upon the request of the client or for any samples which are known to be involved in a legal dispute. As with all analytical data, it is extremely important that this documentation is filled out completely and accurately with every sample transfer. Everyone who handles the COC has the responsibility to check for documentation compliance to the point of their acquisition. If changes need to be made to the form, they shall be made in accordance to the error correction procedure addressed in SOP-QA-109, "Laboratory Notebooks and Documentation." It will be the responsibility of the person who made an error in documentation to correct the error.
- 2. If requested by the client, the COC documentation will begin with the preparation of sampling containers. A form (Figure 1, attached) will be initiated by the person packing the bottle order for shipment to the client. If the delivery of containers is via Lancaster Laboratories Transportation Department, the driver shall sign the form when they relinquish the bottles to the client. Drivers must also sign COC forms when they pick up samples for analysis.
- 3. When samples arrive at the laboratory for analysis, a member of the Sample Administration Group will receive them and sign the external COC form that accompanies the samples, if provided. If the samples were picked up by our Transportation Department, the driver must sign the COC to relinquish the samples to sample administration.
- 4. The Sample Administration Group will track the custody of samples between receipt and entry into the Sample Management System on the SA Receipt Documentation Log (Figure 2, attached). The client's sample designation will be used for identification purposes until a unique Lancaster Laboratories' number is assigned.

Section No. 7 Revision No. 4 Date: 05/12/97 Page 13 of 26

SOP-QA-104.01

Supersedes Date: 12/01/95
Effective Date: MAR 1 4 1997
Page 3 of 16

5. Samples will be entered into the Sample Management System as described in SOP-QA-102, "Sample Log-in." Sample Administration will enter an analysis number for "Laboratory Chain of Custody" if requested. A lab note will print to inform analysts of the need for COC documentation. This note will also be automatically added to the sample labels.

### B. Creating the Internal Chain of Custody

- 1. Sample Administration personnel shall initiate an internal Laboratory Chain of Custody form at the time of sample entry (Figure 3, attached) for each type of container in the sample group. A master list of all chains created will also be initiated for each sample group at the time of entry (Figure 4, attached). The samples will then be relinquished to a sample custodian who will store the samples in an assigned secure location. This change of custody from sample entry to storage shall be documented on the chain, as well as any interim exchanges for rush analysis, preservation, homogenization, or temporary storage in the SA HOLD. The internal COC forms will then accompany the samples from storage to the laboratory for analysis.
- If samples need to be checked out from the Sample Administration Group before Lancaster Laboratories' numbers have been assigned to them, SA will be responsible for starting a COC form. They will note the available header information, the samples being relinquished (documented by the client sample designation), and the reason for transfer.
- 3. After sample entry, the original copy of the external client COC/analysis request form will be filed with Accounts Receivable, to be returned to the client with their invoice. Other copies of the external form will stay within SA to be filed within the client's paperwork file.

Section No. 7 Revision No. 4 Date: 05/12/97 Page 14 of 26

SOP-QA-104.01

Supersedes Date: 12/01/95 Effective Date: MAR 1 4 1997

Page 4 of 16

# C. Documentation of Custody Changes

 An example of how to document changes in sample custody is shown in Figures 3 and 5. Each change of sample custody must be accurately documented in a consistent format. All signatures documenting changes of custody will use the following format:

Signatures: first initial, full last name, employee number

Date: Month/day/year

Time: Documented as military time

Ink: Black ink is preferred, red ink and pencil are not acceptable

a. When sample support releases samples to an analyst they must:

Note the sample number(s) released, and sign the released by column of the chain.

b. When an analyst receives samples from sample support they must:

Sign the received by column, note the date and time samples are received and note the reason why they are taking the samples (reason for change of custody).

c. When an analyst returns samples to sample support they must:

Note all sample numbers being returned, sign the released by column, and note time and date of return.

d. When sample support receives samples from an analyst they must:

Sign the received by column and note the reason for sample transfer.

Section No. 7 Revision No. 4 Date: 05/12/97 Page 15 of 26

SOP-QA-104.01

Supersedes Date: 12/01/95 Effective Date: MAR 1 4 1997

Page 5 of 16

- 2. Sample handling should be kept to a minimum. Analysts requiring use of a sample will requisition it through the computer requisition program. During the hours when sample support is manned by sample custodians, a custodian will receive the computerized requisition and remove the sample from storage. The custodian will ensure that the bottle type listed on the COC form matches the bottle type being distributed. It will be the shared responsibility of the analyst and sample custodian to insure that forms are signed, dated, and reason for sample transfer are recorded with each change of custody, as directed by Item C1 above.
- Each specific test that an analyst performed in conjunction with the
  associated sample number(s) must be accurately documented by the analyst
  before the samples are returned to a sample custodian in the sample storage
  area.
- 4. When an analyst requires the use of samples when a sample custodian will not be on duty, they must requisition samples earlier in the day or on the previous day. These samples and associated COCs will be pulled by a sample custodian and placed in the locked SA HOLD storage area. The sample custodian will note on the COC the change in transfer to the SA HOLD in addition to the time, date, and the sample numbers. The analyst picking up the samples will document the specific samples being checked out, record SA HOLD in the "Released by" column, sign the Received by column, note the time, date and reason for transfer. When the analyst returns the samples to the SA HOLD, they must sign the samples back into the SA HOLD.
- 5. The following changes of custody will be handled in the following manner:
  - a. Documentation is required for all shift changes. Signatures involving transfers from one shift to another shall be the responsibility of the analyst who originally acquired the samples from sample support.

Section No. 7 Revision No. 4 Date: 05/12/97 Page 16 of 26

SOP-QA-104.01 Supersedes Date: 12/01/95 Effective Date: MAR 1 4 1997 Page 6 of 16

- b. Occasionally a sample container will be needed for analysis by an analyst in a department while it is in the custody of an analyst in another department. It will be the responsibility of the first person who received the sample to note on the COC the specific sample numbers requested by the second person, and to sign the released by column. The second person will sign the received by column and note the time, date, and reason for sample transfer. After the second person is finished with the sample, the sample will be returned back to the first person or to the sample storage area.
- c. In situations where a sample group must be split between departments working on different analyses, a supplemental COC must be initiated by the Sample Support Group. The supplemental chain will be used to accompany that portion of the sample group which is needed by a second department, when another department has part of the sample group and the COC for the entire group. This supplemental COC will be created only when absolutely necessary to minimize paperwork and confusion. This chain must also be documented on the master list of chains initiated for the sample group.
- d. Some original samples are released by Sample Support or Sample Administration to be stored in other areas of the laboratory (e.g. GC/MS Volatiles, Foods, Pharmaceuticals). During this time they may be accessed by several people in that area. Each of these people must note the specific sample numbers in their custody in addition to date, time, and reason for removal from storage. An example of a COC is attached as Figure 6.

It will be the responsibility of the department who held the samples to assure that all necessary, signatures, dates, times, and reasons for sample custody are noted on the COC forms. It is also very important to return all samples and COCs to storage as soon as possible after data verification, because the chains may be required for a client data packages.

Section No. 7 Revision No. 4 Date: 05/12/97 Page 17 of 26

SOP-QA-104.01

Supersedes Date: 12/01/95 Effective Date: MAR 1 4 1997

Page 7 of 16

e. If COC samples are stored in other areas of the laboratory or in a specific department, they must be stored in a locked area. When samples are taken from a departmental storage area, the released by column of the COC is documented as "department XX storage." If samples are returned to this area when complete the received by column will be noted as department XX storage.

# D. Additional Chain-of-Custody Issues

- Analysts in possession of samples shall remove the aliquot required for their analysis and return the samples to the Sample Support Group with a minimum of delay. During this time of possession, samples must fall under the definition of sample custody.
- 2. If additional containers of the sample are created (e.g. subsamples, extracts, distillates, leachates, digests, etc.), an additional COC form must be created by the department if they do not document this information on the original COC form (Figure 5, attached). This form will be marked with the container type and will be initiated to accompany the new sample container. Each department in the lab has specifically designed COC forms which will be used if new containers are created. All changes of custody involving handling of new containers in the department (e.g. analysis, storage, vials on instruments, etc.) will be documented on the departmental specific COC form or on the original COC form. Any specific handling or documentation requirements for departmental chains can be described in a departmental SOP.

#### E. When Sample Analysis is Complete

After sample analysis, samples shall be returned to the Sample Support
Group as soon as possible. Original COC forms shall also be returned with
the samples and this change of custody noted. At this time it will be the

Page 18 of 26

SOP-QA-104.01

Supersedes Date: 12/01/95 Effective Date: MAR 1 4 1997

Page 8 of 16

responsibility of the Sample Support Group to review the COC forms to ensure that all documentation on the forms is complete before they file the forms in their area. Sample custodians will not return a sample to its assigned storage location without signing the accompanying chain and performing this completeness check. All chains should either end with a note of "All Sample Consumed," "Discard" or "Storage" for the final reason of transfer.

- 2. All completed COC forms for the original sample containers will be retained in files within Sample Support. The Data Package Group will retrieve these forms so a copy can be included in the data package. All departmental created COC forms will be collected by the department's data package group so a copy can be included in the data package. These forms will not be returned to the Sample Support Group since these sample containers will not be returned to the Sample Support Group. The original copy of all COC forms will be retained on file by the laboratory.
- 3. All personnel who handle sample containers shall make every attempt to ensure that all changes of custody are accurately and completely documented. Disciplinary action may be taken for employees who fail to comply with these important requirements.
- 4. In the event that a signature or other information is inadvertently not recorded on a COC form, the Sample Support and Data Package Groups in conjunction with the technical centers shall determine what information is missing by checking computer requisition records, raw data, or the sample support work schedule. The responsible party shall add the missing information or make the necessary correction at the bottom of the COC form, in addition to noting the situation that caused the error in documentation. The person making this note needs to sign and date the information using the current date. Any errors in COC documentation that cause noncompliances must be noted in the case narrative of the sample data package. Examples of specific cases are on file in the data package department.

Section No. 7 Revision No. 4 Date: 05/12/97 Page 19 of 26

SOP-QA-104.01

Supersedes Date: 12/01/95 Effective Date: MAR 1 4 1997

Page 9 of 16

# Personnel Training & Qualification:

Training for this procedure consists of reading SOP-QA-104. Supervisory review of all COC documentation should be done until the trainer is satisfied that proficiency has been achieved. Training of all laboratory personnel is the responsibility of the group leader. Documentation that this training has been completed must be kept in the training records.

# **Revision Log:**

Initiated Date: 03/87

Ver. # Effective Date Change
00 12/01/95 Previous Issue
01 MAR 1 4 1997 Major changes are as follows:

- Training section added.
- Examples of SA Receipt Documentation Log and Metals Locked Storage COC updated.
- Section E.1., Option to end chain with "All Sample Consumed" added.

SOPQA104.DOC 021997

Prepared by:

972 Date:

2/28/97

Approved by:

Data:

3/5/97

Approved by:

Date:

3|5|97

For LLI use only

Lancaste	r Laboratories													^	icct.			_ Sample #	•			-
<u> </u>		lease pri	nt. Instru	ctio	ns or	rev	erse :	side	com	espon	d with			nbers.								
Clent:		Acct. #:_				_	trix (	<u> </u>	1		ઉ	$\sum$		A	nalys	es Re	quest	ed /	For I	LLI use	-	7
Project Name/#:		PWSID #:				<u>.</u>				ĺ	/	/ /	′,	/ /	/ ,	/ /	/ /	/	SCR 0:_			
Project Manager:		P.O.#					Oect	50	·												١,	(6
Sampler;		Quote #:				• .	ង្គី ស		Containers			/	/					/ /			19043	
							C Potable	3.4	70	/	//		′ /	/ /	/ /	/ /	/ /				183	ı i
Name of state where samples w  Sample Identification		Date Collected	Time Collected	de S	Comp	3	Water	Other	Total #								_	Remarks		 	Temperature Temperature	
				L																		]
				L	Ш			<u> </u>	L						<u> </u>	<u> </u>	L				_ _	1
				L	Ш				_		إـــا		<u> </u>	↓	<del> </del>	<u> </u>	ļ	ļ		·		1
·		*****		L	Ш			<u> </u>	_	<u> </u>			ļ	<b>↓</b>	<u> </u>	. <b>.</b>	↓					]
				ļ.,	Ш				ļ				<u> </u>	1	<del> </del>	↓	ļ					4
			ļ · ·	<u> </u>	H				<u> </u>	<u> </u>	ļ		ļ	<b></b>	╄	ļ	ـــــ	<del> </del>				4
			<u> </u>	-	Н				_	<u> </u>			<u> </u>	╄	<del> </del>	1_	<u> </u>	<u> </u>				-
· · · · · · · · · · · · · · · · · · ·			ļ	ļ	Н			<u> </u>	<u> </u>	-				1-	╂	<del> </del>	├-	<del> </del>				4
<del></del>			<del> </del>	-	Н	-		<del> </del>	_				_	-	+-	┨	-	<del> </del>				4
Turnaround Time Requeste (Rush TAT is subject to LLI approval as		ormal R	lush	Ret	inquis	hedi	<b>)</b> 3 <b>y</b> :	L	<u> </u>	<u>l</u>		T 64	ate	Time	Recei	ved by:	ــــــــــــــــــــــــــــــــــــــ	<u> </u>		Date	Time	ţ
Date results are needed: Rush results requested by (please Phone #:				Ret	nquis	hed t	y:					Di	ote	Time	Recei	ved by:			$\neg$	Date	Time	1
Data Package Options when			Complete?	Ret	nquis	hed I	y:					D	ile	Time	Recei	ved by:				Date	Time	7
Type I (Tier I) GLP	Site-specific QC required If yes, indicate QC sample and	7 Yes		Ret	nqus	hed I	y.					D	ne	Time	Recei	ved by:				Date	Time	
Time M (NII Day) Dall	nternal Chain of Custod	·-			nquis	hed I	γ.					D,	rie	Time	Recei	ved by:				Date	Time	٦

2425 New Holland Price, PD Box 12425, Lancester, PA 17605-2425 (717) 656-2300 Copies: White and yellow should accompany samples to Lancester Laboratories. The pink copy should be retained by the Chent

SOP-QA-104.01
Supersedes Date: 12/01/95
Effective Date: MAR 1 4 199
Page 10 of 16

Section No. 7
Revision No. 4 Date: 05/12/97

Page 20 of 26

#### **DIRECTIONS FOR COMPLETING THIS FORM**

(1) Client: Your company's name

Acct. #: Your account number with Lancaster Laboratories

Project Name/#: The way your company refers to the work involved with these samples. You may want to include project location as part of the description.

PWSID: Potable Water Source ID#

Project Manager: The person at your company responsible for overseeing the project

P.O. #: Your company's purchase order number

Sampler: The name of the person who collected the samples

Quote #: The reference number that appears on your quote (if Lancaster Laboratories gave you a number)

State where sample was collected: Please indicate where the sample was taken, e.g., Pa., N.J., etc.

(2) Sample Identification: The unique sample description you want to appear on the analytical report

Date Collected/Time Collected: When the sample was collected

- (3) Grab: Check here if sample was taken at one time from a single spot. Composite: Check here if samples were taken from more than one spot, or periodically, and combined to make one sample.
- (4) Matrix: Check the type of sample you are submitting. If it is a water sample, please indicate if it is a potable water or if it is an NPDES

Number of Containers: Indicate the total number of containers for each sampling point.

- (5) Analyses Requested: Write the name of each analysis for an abbreviation of it) here, and use the catalog number that appears at the beginning of each line in the Schedule of Services. Be sure to indicate which analyses are to be performed on which samples.
- (6) Remarks: List special instructions about the sample here (e.g., hazardous elements, high levels of analyte, etc.). The space can also be used (if needed) for listing additional analyses.
- (7) Turnaround time Requested: Circle Normal if you want routine TAT, which is usually within 10-15 days. If you need your results faster, call ahead to schedule Rush work.

Rush Results Requested by: Circle Fax or Phone and include the number.

(8) Data Package Options: Call our Client Services Group (717-656-2301) if you have questions about these choices.

SDG Complete? Indicate Yes if this is a complete sample delivery group or No if you will be submitting additional samples to be included in the same data package.

Note: We need to have one quality control (QC) sample for every 20 samples you send, if you are requesting site-specific OC. Please give us this sample in triplicate volume and identify it by writing "QC" in the Remarks column.

The internal chain of custody is a hand-to-hand documentation recording a sample's movement throughout the company. We routinely start a chain of custody for data-package samples unless we are told otherwise. There is a \$25 per sample charge for the chain-of-custody documentation.

(9) Relinquished by/Received by: The form must be signed each time the sample changes hands. We can supply chain-of-custody seals for the outside of your packages if you require them.

Thank you for using Lancaster Laboratories.

Please call our Client Services Group (717-656-2301) if you have any questions about completing this form.

Figure 1 - Continued

SOP-un Supersedes Date Effective Date: Page 11 of 16

MAR 14

upersedes Date OP-QA-104.01 Page Date: 21 of 26

Section No. 7 Revision No.

Section No. 7 Revision No. 4 Date: 05/12/97 Page 22 of 26

SOP-QA-104.01

Supersedes Date: 12/01/95 Effective Date: MAR 1 4 1997

Page 12 of 16

# Figure 2



# Sample Administration Receipt Documentation Log

Client/Project: XVZ Asso	ciates / well-torin	p COC Seal: Pr	esent Not Present on cooler					
Date of Receipt: 11/27		,	roken / htact					
Time of Receipt: 1350	•	Package Chilled Not Chilled						
Source Code: 60		Unpacker Emp. No.: 2/0						
	Temperature (	of Samples						
#1			#2					
Thermometer ID: 123		Thermometer I	D:					
Corrected Temp.: NA		Corrected Tem	р.:					
Temp. Bottle / Surface Temp	o.	Temp. Bottle	Surface Temp.					
Wet Ice / Ice Pa	cks	Wat Ice / Dry	Ice / Ice Packs					
Ice Present? (Y) N		Ice Present?	Y / N					
#3			#4					
Thermometer ID:		Thermometer I	D:					
Corrected Temp.:		Corrected Tem	p.;					
Temp. Bottle / Surface Tem	p.	Temp. Bottle	Surface Temp.					
Wet Ice / Dry Ice / Ice Pac	ks .	Wet Ice / Dry	ice / ice Packs					
Ice Present? Y / N		Ice Present?	Y / N					
Paperwork Discrepancy/Unpag	king Problems: Book	ten 40 .	ne vial, chient 1D					
123-01 Ch	int called	1610	11/22/05					
165-01	ec. Chius	1010	1421/1)					
			<del></del>					
	Sample Administration	n Chain of C	ustody					
Name	Date	Time	Reason for Transfer					
L'ant.	11/27/95	1600	Unpacking					
A Hulchism	11/27/25	1615 (	Place in Storage or Entry					
D. Nestural	11/22/55	1800	Remove from Storage Entry					
	<del></del>		Place in Storage or Entry					
			Entry					

Section No. 7 Revision No. 4 Date: 05/12/97 Page 23 of 26

SOP-QA-104.01

Supersedes Date: 12/01/95 Effective Date: MAR 1 4 1997

Page 13 of 16



Client/Project: XYZ Assisted

# Figure 3 Locked Storage Chain of Custody Original Sample

Preservative:	<u> C </u>		Matrix:	Votes)_	sdg: <u>XYZC</u>	)/
Sample # Range of E	ntry Group:	2420	632-39	<b>)</b>	Bottle Type: ##	+38 L
Sample Number(s) in Custody	Released By	Received By	Date of Transfer	Time of Transfer	Reason for Change of Custody	Dist., Extr., or Digest Chain Created (X)
<i>24</i> 20638-39	D. 208 Meslund	55 Looge	11/27/95	· -	Entry & Storage	
2420638-39	55 Storage	B. 705 Weaver	11/28/95	700	Hemole from 55 Storage	
242 0638-39	B. 705 Weaven	dipt 21 Strage	11/28/95	715	VOA Storage	
2420638-39	dept 21 Storage	K. 3% Witman	11/29/95	1315	VOA analysis	X
2420638-39	K. Witman 396	d. 53 Toylor	11/29/95	1700	VOA aralyst Shift Charge	
2420638-39	f. 513 Taylor	dept 21 Storage	11/29/95	2100	VOA Storage	
a420638:39	dept 21 Storage	C. 24.	12/3/95	<i>80</i> 0	Transfer to 55 Storage	
2420638-39	C. Oypro	2. 630 Level V	12/3/95	815	Storage	
·						

_		
<sup>D</sup> age	~ f	
auc	of	

Section No. 7 Revision No. 4 Date: 05/12/97 Page 24 of 26

SOP-QA-104.01

Supersedes Date: 12/01/95 Effective Date: Page 14 of 16 MAR 1 4 1997

<b>4</b> 12	Lancaster	Laboratories
717	3425 New Helland Pile	· Lancaster, PA 17601-5994

# Figure 4 Master List of Chain of Custodies

Client/Project: XYZ amociates	·	* · · · · · · · · · · · · · · · · · · ·
Sample # Range of Entry Group: 24200	3 <i>3 -3</i> 9	
SDG: XYZ01 Matrix: Liq	uid Solid Mixed Othe	er
Original Sa	mple Chains	
Bottle Type	Started By	Date Started
40 ml bloss Vial (#38)	D. Meclind 208	11/27/95
1000 ml amber Hinser (*45)		
1000 ml Plastic (#09)		
1000 ml amber Glass (*29)		<u> </u>
	<u> </u>	
Supplemen	ntal Chains	
Bottle Type	Started By	Date Started
77	C. Opa 266	11/27/95
ત્રા	C. Eyard 266	11/27/95
	U	
Extraction, Digesti	on, Distillates, Etc.	
Bottle Type	Started By	Date Started

Section No. 7 Revision No. 4 Date: 05/12/97 Page 25 of 26

SOP-QA-104.01

Supersedes Date: 12/01/95 Effective Date: MAR 1 4 1997

SDG:

XYZOI

Trial No: 2 (If not 1, fill in)

Page 15 of 16

# Figure 5



Client/Project: XYZ Associates

(Metals )

Sample #: 2420632-9

Digest Type (circle one): Hg

## Locked Storage Chain of Custody Metals

Hydrides

Batch No:	9530	5 18	49	0	04
Sample Number(s) in Custody	Released By	Received By	Date of Transfer	Time of Transfer	Reason for Change of Custody
2420632,4,6	S. Correa/523  J. Granett/428  ICl Stunge	J. Garrett /428	12/1/95	163)	Mital pres Shiff change TCP Storage TCP Strage.
2420632,4,6	J. Garrett/425	ICPStorage	12/1/95	1920	ICP Storage
2420632,4,6	ICIStonge		12/1/95	2115	Il Analysis
2420632,46	D. Sachett/	ICP Stray	12/1/95	2/35	ICP Shaye.

Section No. 7 Revision No. 4 Date: 05/12/97 Page 26 of 26

SOP-QA-104.01

Supersedes Date: 12/01/95 Effective Date: MAR 1 4 1997

Page 16 of 16



Client/Project: <u>Smith Pharmaceutical</u>

# Figure 6

# Pharmaceutical Locked Storage Chain of Custody Original Sample

Preservative:	N/A	Matrix: <u>//</u>	weis		<del></del>
Sample # Range of E					<del></del>
Bottle Type:	40 ml Vial				
Sample Number(s) in Custody	Released By	Received By	Date of Transfer	Time of Transfer	Reason for Change of Custody
a420320 30	S. 014 Carrethers	m. 589 Coho	12/1/95	1300	Entry/Transfer to Thaim. Storage
a4 <i>a</i> 03a0-30	m. 589 Coho	Storage	18/1/95	1315	0
2420330-30	Storage	e. 5721	12/3/95	800	Storage PN analysia
3420320-30	etc. 572	Storage	12/3/95	1000	Gorage
A4A03A0 :30	Storage	D. 330 Wright	12/5/95	930	GC assayl analysis Storage
a4a03a0-30	D. 330 Wright	Storage	12/5/95	1400	Storage

Section No. 8 Revision No. 2 Date: 07/25/95 Page 1 of 5

# 8. Calibration Procedures

Procedures for initial calibration and continuing calibration verification are in place for all instruments within the laboratory. The calibrations generally involve checking instrument response to standards for each target compound to be analyzed. The source and accuracy of standards used for this purpose are integral to obtaining the best quality data. Standards used at Lancaster Laboratories are purchased from commercial supply houses either as neat compounds or as solutions with certified concentrations. The accuracy and quality of these purchased standards is verified through documentation provided by these commercial sources. Most solutions and all neat materials require subsequent dilution to an appropriate working range. All dilutions performed are documented and the resulting solution is checked by obtaining the instrument response of the new solution and comparing with the response to the solution currently in use. Any discrepancies between the responses are investigated and resolved before the new solution is used. Each standard is assigned a code which allows traceability to the original components. The standard container is marked with the code, name of solution, concentration, date prepared, expiration date, and the initials of the preparer. Shelf life and storage conditions for standards are included in the standard operating procedures and old standards are replaced before their expiration date.

Each instrument is calibrated with a given frequency using one or more concentrations of the standard solution. As analysis proceeds, the calibration is checked for any unacceptable change in instrument response. If the calibration check verifies the initial response, the analysis proceeds. If the calibration check indicates that a significant change in instrument response has occurred, then a new calibration is initiated. If necessary, maintenance may be performed prior to the recalibration.

Calibration records are usually kept in the form of raw data with the other instrument printouts. In cases where no data system is used, calibration data is manually recorded in notebooks. Any maintenance or repair is also recorded in a notebook. The information recorded either in the notebooks or on the instrument printout includes the date, instrument ID, employee name and/or identification number, and concentration or code number of standard.

Section No. 8 Revision No. 2 Date: 07/25/95 Page 2 of 5

The frequency of calibration and calibration verification, number of concentrations used, and acceptance criteria for each of the instruments to be used are listed on Table 8-1. In addition to checking the instrument response to target compounds, the GC/MS units are checked to ensure that standard mass spectral abundance criteria are met. Prior to each calibration, instruments being used for volatile compound analysis are tuned using bromofluorobenzene (BFB) and instruments being used for semivolatile analysis are tuned using decafluorotriphenylphosphine (DFTPP). The key ions and their abundance criteria are listed in Table 8-2.

Section No. 8 Revision No. 2 Date: 07/25/95 Page 3 of 5

			Table 8-1			
			alibration	Continu		ration Verification
Instrument	Frequency	# Std Conc	Acceptance Criteria	Frequency	# Std Conc	Acceptance Criteria
GC/MS Volatiles	After C-cal fails	5	Specified compounds must meet contract minimum RRF criteria and max %RSD of <20.5%	Every 12 hours	1	Specified compounds must meet contract minimum RRF criteria and max %D of ≤25%
GC/MS Semivolatiles	After C-cal fails	5	Specified compounds must meet contract minimum RRF criteria and max %RSD of ≤20.5%	Every 12 hours	1	Specified compounds must meet contract minimum RRF criteria max %D of ≤25%
GC Pesticides	After C-cal fails	3	%RSD for compounds ≤20% (alpha-BHC and delta-BHC ≤25%) %RSD for surrogates ≤30%	Every 12 hours	1	INDA&B/PEM alternate every 12 hours with %D ≤25. Degradation for DDT, endrin ≤20%, combined ≤30%
Flame AA	Each new run	5	Independent calibration within ±10%	Every 10 samples	1	Same as initial
Cold Vapor AA	Each new run	5	Independent calibration verification within ±20%	Every 10 samples	1	Same as initial
ICP	Each new run Max. 60 samples- run	2	Independent calibration verification within ±10%	Every 10 samples	1	Same as initial
Graphite Furnace AA	Every new run	5	Independent calibration verification within ±10%	Every 10 samples	1	Same as initial
Autoanalyzer (cyanide)	Daily	5	Correlation coefficient >0.995	Every 10 samples	1	±10% of original response
Balance	Daily	4	±.5%	N/A	N/A	N/A

# **Abbreviations**

- RRF Relative response factor
  - %RSD Percent relative standard deviation
  - %D Percent difference
  - RPD Relative percent difference
  - C-cal Continuing calibration
- Flame AA Flame atomic absorption spectrophotometer
  - ICP Inductively coupled plasma spectrophotometer
- Graphite Furnace AA Graphite furnace atomic absorption spectrophotometer

Section No. 8 Revision No. 2 Date: 07/25/95 Page 4 of 5

For volatiles, up to two compounds may be outside criteria providing the RRF is ≥0.010 and %RSD ≤40%.

For semivolatiles, up to four compounds may be outside criteria providing the RRF is ≥0.010 and %RSD ≤40%.

For both volatile and semivolatile compounds with no established RRF criteria, the minimum RRF is ≥0.010.

For pesticides, up to two target compounds may have %RSD >20% but ≤ 30%.

Section No. 8 Revision No. 2 Date: 07/25/95

Page 5 of 5

Table 8-2							
Mass	Ion Abundance Criteria						
BFB Key Ion Abundance Criteria:							
50	8% to 40% of mass 95						
75	30% to 66% of mass 95						
95	base peak, 100% relative abundance						
96	5% to 9% of mass 95						
173	less than 2% of mass 174						
174	50% to 120% of mass 95						
175	4% to 9% of mass 174						
176	93% to 101% of mass 174						
177	5% to 9% of mass 176						
DFTPP Key lons	and Ion Abundance Criteria:						
51	30% to 80% of mass 198						
68	less than 2% of mass 69						
69	mass 69 relative abundance						
70	less than 2% of mass 69						
127	25% to 75% of mass 198						
197	less than 1% of mass 198						
198	Base peak, 100% relative abundance						
199	5% to 9% of mass 198						
275	10% to 30% of mass 198						
365	greater than .75 of mass 198						
441	Present but less than mass 443						
442	40% to 110% of mass 198						
443	15% to 24% of mass 442						

Page 1 of 9

# 9. Analytical Procedures

The analytical procedures to be used for organics are those described in the USEPA CLP Organics SOW OLM03.2 or most recent version and are designed to analyze water, sediment, and soil for the organic compounds on the Target Compound List (TCL). The inorganics procedures are those indicated in the USEPA CLP Inorganics SOW ILM04.0 for the preparation and analysis of water, sediment, and soil samples for the elements on the Target Analyte List (TAL). Copies of the analytical procedures are located in the laboratory and are available for use by analysts. Copies of analytical methods are available upon request.

<u>Volatiles</u> - This method determines the concentration of TCL volatile (purgeable) organics. The analysis is based on purging the volatiles onto a Tenax/silica gel trap, desorbing the volatiles onto a gas chromatographic column which separates them and identifying the separated components with a mass spectrometer. (GC/MS Method.)

<u>Semivolatiles</u> - This method determines the concentration of semivolatile organic compounds that are separated into an organic solvent and are amenable to gas chromatography. The method involves solvent extraction of the sample to isolate analytes and GC/MS analysis to determine semivolatile (BNA) compounds present in the sample.

<u>Pesticides</u> - This method determines the concentration of TCL organochloride pesticides and polychlorinated biphenyls. The procedure includes solvent extraction of the sample, analysis of the extract on a gas chromatograph/electron capture detector (GC/EC) using a megabore capillary column, and confirmation on a GC/EC using a second megabore capillary column. If the compound concentration is sufficient, confirmation may be done by GC/MS upon request.

Inductively Coupled Plasma (ICP) - This is a technique for the simultaneous determination of elements in solution after acid digestion. The basis of the method is the measurement of atomic emission by an optical spectroscopic technique. Characteristic atomic line emission spectra are produced by excitation of the sample in a radio frequency inductively coupled plasma. Because of the high temperature of the plasma, it is especially useful for refractory metals.

Section No. 9 Revision No. 5 Date: 05/12/97 Page 2 of 9

The Trace ICP is the same technique as the ICP listed above except for the orientation of the plasma (horizontal instead of vertical) and upgraded optical and sample introduction systems, resulting in instrument detection limits approximately a magnitude lower than the traditional ICP.

Graphite Furnace Atomic Absorption (GFAA) - This is a method of analysis designed to detect trace amounts of the analyte through electrothermal atomization. Samples are digested before analysis. The graphite furnace is an AA spectrophotometer that heats the sample within a graphite tube using an electrical current (i.e., flameless furnace) and measures the absorption of specific metallic elements at discrete wavelengths.

<u>Flame Atomic Absorption</u> - This method is also suited to metals analysis. A solution of the sample to be analyzed is sprayed into a flame which generates sufficient heat to decompose the sample into its constituent atoms directly in the optical path. The difference in light intensity is measured at specific wavelengths using a spectrophotometer.

<u>Cold Vapor Atomic Absorption</u> - Organic mercury compounds are oxidized and the mercury is reduced to the elemental state and aerated from solution in a closed system. The mercury vapor passes through a cell positioned in the light path of an AA spectrophotometer and absorbance (peak height) is measured.

<u>Total Cyanide Analysis</u> - Cyanide, as hydrocyanic acid, is released from cyanide complexes by means of a reflux-distillation operation and absorbed in a scrubber containing sodium hydroxide solution. The cyanide ion in the absorbing solution is then determined colorimetrically.

<u>Percent Moisture</u> - A known sample weight is placed in a drying oven maintained at 103° to 105°C for 12 to 24 hours. The sample is reweighed after drying and this value is divided by the original weight. The result is used to calculate analytical concentration on a dry-weight basis.

Page 3 of 9

Target Compound List (TCL) and Contract Required Quantitation Limits (CRQL)*									
		Water	· (μg/L)	Low S	oil (μg/kg)				
Volatiles	CAS Number	Quant. Limit	J-Value	Quant. Limt**	J-Value**	Med Soil <sup>a</sup> (μg/kg)**			
Chlorornethane	74-87-3	10	3	10_	2	1200			
Bromornethane	74-83-9	10	3	10	3	1200			
Vinyl Chloride	75-01-4	10	2	10	2	1200			
Chloroethane	75-00-3	10	3	10	3	1200			
Methylene Chloride	75-09-2	10	2	10	2	1200			
Acetone	67-64-1	10	6	10	7	1200			
Carbon Disulfide	75-15-0	10	3	10	3	1200			
1,1-Dichloroethene	75-35-4	10	1	10	2	1200			
1,1-Dichloroethane	75-34-3	10	2	10	1	1200			
1,2-Dichloroethene (total)	540-59-0	10	2	10	2	1200			
Chloroform	67-66-3	10	1	10	11	1200			
1,2-Dichloroethane	107-06-02	10	2	10	2	1200			
2-Butanone	78-93-3	10	3	10	7	1200			
1,1,1-Trichloroethane	71-55-6	10	11	10	1	1200			
Carbon Tetrachloride	56-23-5	10	1	10	11	1200			
Bromodichloromethane	75-27-4	10	1	10	2	1200			
1.2-Dichloropropane	78-87-5	10	11	10	3	1200			
cis-1,3-Dichloropropene	10061-01-5	10	1	10	1	1200			
Trichloroethene	79-01-6	10	1	10	1	1200			
Dibromochloromethane	124-48-1	10	2	10	1	1200			
1,1,2-Trichloroethane	79-00-5	10	_2	10	2	1200			
Benzene	71-43-2	10	1	10	1	1200			
trans-1.3-Dichloropropene	10061-02-6	10	1	10	1	1200			
Bromoform	75-25-2	10	1	10	1	1200			
4-Methyl-2-pentanone	108-10-1	10	5	10	3	1200			
2-Hexanone	591-78-6	10	7	10	3	1200			
Tetrachloroethene	127-18-4	10	1	10	11	1200			
Toluene	108-88-3	10	2	10	1	1200			
1 1,2,2-Tetrachloroethane	79-34-5	10	2	10	1	1200			

Page 4 of 9

Target Compound List (TCL) and Contract Required Quantitation Limits (CRQL)\* Water (µg/L) Low Soil (µg/kg) Med Quant. Soil<sup>a</sup> Quant. **Volatiles CAS Number** Limit J-Value Limt\*\* J-Value\*\*  $(\mu g/kg)^{**}$ Chlorobenzene 108-90-7 10 1 1200 10 1 2 1 Ethyl Benzene 100-41-4 10 10 1200 10 1 1200 Styrene 100-42-5 10 1 1 1 1330-20-7 10 10 1200 Xylene (total)

J-values are evaluated annually and are subject to change.

<sup>\*</sup>Specific quantitation limits are highly matrix dependent. The quantitation limits listed herein are provided for guidance and may not always be achievable.

<sup>\*\*</sup>Quantitation limits and J-values listed for soil/sediment are based on wet weight. The quantitation limits calculated by the laboratory for soil/sediment, calculated on dry weight bases as required by the contract, will be higher.

<sup>&</sup>lt;sup>a</sup>The J-value for the medium-level soil analysis can be determine by multiplying the low-level soil J-value by a factor of 125 and then rounding according to CLP protocol.

Page 5 of 9

Target Compound List (TCL) and Contract Required Quantitation Limits (CRQL)*									
		Water	r (μg/L)	Low S					
Semivolatiles	CAS Number	Quant. Limit	J-Value	Quant. Limt**	J-Value**	Med Soil <sup>a</sup> (μg/kg)**			
Phenol	108-95-2	10	1	330	67	10000			
bis(2-Chloroethyl)ether	111-44-4	10	1	330	67	10000			
2-Chlorophenol	95-57-8	10	1	330	67	10000			
1,3-Dichlorobenzene	541-73-1	10	1	330	33	10000			
1,4-Dichlorobenzene	106-46-7	10	1	330	33	10000			
1,2-Dichlorobenzene	95-50-1	10	1	330	33	10000			
2-Methylphenol	95-48-7	10	2	330	67	10000			
2,2'-oxybis(1-Chloropropane)	108-60-1	10	1	330	33	10000			
4-Methylphenol	106-44-5	10	2	330	67	10000			
N-Nitroso-di-n-dipropylamine	621-64-7	10	1	330	33	10000			
Hexachloroethane	67-72-1	10	1	330	33	10000			
Nitrobenzene	98-95-3	10	1	330	33	10000			
Isophorone	78-59-1	10	1	330	33	10000			
2-Nitrophenol	88-75-5	10	1	330	33	10000			
2,4-Dimethylphenol	105-67-9	10	2	330	33	10000			
bis(2-Chloroethoxy)methane	111-91-1	10	1	330	33	10000			
2,4-Dichlorophenol	120-83-2	10	1	330	33	10000			
1,2,4-Trichlorobenzene	120-82-1	10	1	330	33	10000			
Naphthalene	91-20-3	10	1	330	33	10000			
4-Chloroaniline	106-47-8	10	1	330	33	10000			
Hexachlorobutadiene	87-68-3	10	1	330	33	10000			
4-Chloro-3-methylphenol	59-50-7	10	1	330	33	10000			
2-Methylnaphthalene	91-57-6	10	1	330	67	10000			
Hexachlorocyclopentadiene	77-47-4	10	1	330	100	10000			
2,4,6-Trichlorophenol	88-06-2	10	1	330	33	10000			
2,4,5-Trichlorophenol	95-95-4	25	1	830	33	25000			
2-Chloronaphthalene	91-58-7	10	1	330	33	10000			
2-Nitroaniline	88-74-4	25	1	830	33	25000			

Section No. 9 Revision No. 5 Date: 05/12/97 Page 6 of 9

Targe	et Compound List (TCL) and
Contract Req	quired Quantitation Limits (CRQL)*

Contract Required Quantitation Limits (CRQL)*								
		Water (μg/L)		Low S	oil (μg/kg)			
Semivolatiles	CAS Number	Quant. Limit	J-Value	Quant. Limt**	J-Value**	Med Soil <sup>a</sup> (μg/kg)*'		
Dimethylphthalate	131-11-3	10	1	330	33	10000		
Acenaphthylene	208-96-8	10	1	330	33	10000		
2,6-Dinitrotoluene	606-20-2	10	1	330	33	10000		
3-Nitroaniline	99-09-2	25	1	830	67	25000		
Acenaphthene	83-32-9	10	1	330	33	10000		
2,4-Dinitrophenol	51-28-5	25	11	830	33	25000		
4-Nitrophenol	100-02-7	25	1	830	67	25000		
Dibenzofuran	132-64-9	10	1	330	33	10000		
2,4-Dinitrotoluene	121-14-2	10	1	330	33	10000		
Diethylphthalate	84-66-2	10	1	330	33	10000		
4-Chlorophenyl-phenyl ether	7005-72-3	10	11	330	33	10000		
Fluorene	86-73-7	10	1	330	33	10000		
4-Nitroaniline	100-01-6	25	2	830	33	25000		
4,6-Dinitro-2-methylphenol	534-52-1	25	1	830	33	25000		
N-nitrosodiphenylamine	86-30-6	10	1	330	33	10000		
4-Bromophenyl-phenylether	101-55-3	10	1	330	33	10000		
Hexachlorobenzene	118-74-1	10	_11	330	33	10000		
Pentachlorophenol	87-86-5	25	2	830	100	25000		
Phenanthrene	85-01-8	10	1	330	33	10000		
Anthracene	120-12-7	10	1	330	33	10000		
Carbazole	86-74-8	10	1	330	33	10000		
Di-n-butylphthalate	84-74-2	10	1	330	100	10000		
Fluoranthene	206-44-0	10	1	330	33	10000		
Pyrene	129-00-0	10	1	330	33	10000		
Butylbenzylphthalate	85-68-7	10	1	330	33	10000		
3,3'-Dichlorobenzidine	91-94-1	10	4	330	33	10000		
Benzo(a)anthracene	56-55-1	10	1	330	33	10000		
Chrysene	281-01-9	10	1	330	100	10000		
bis(2-Ethylhexyl)phthalate	117-81-7	10	1	330	33	10000		

Section No. 9 Revision No. 5 Date: 05/12/97 Page 7 of 9

Target Compound List (TCL) and Contract Required Quantitation Limits (CRQL)\* Water (µg/L) Low Soil (µg/kg) Med Soil Quant. Quant. **Semivolatiles CAS Number** J-Value Limt\*\* J-Value\*\* Limit (μg/kg)\*\* Di-n-octylphthalate 117-84-0 10 330 10000 1 33 10 1 10000 Benzo(b)fluoranthene 205-99-2 330 33 207-08-9 10 1 330 33 10000 Benzo(k)fluoranthene 10 1 330 33 10000 50-32-8 Benzo(a)pyrene 193-39-5 1 10000 Indeno(1,2,3-cd)pyrene 10 330 33 Dibenz(a,h)anthracene 53-70-3 10 1 330 33 10000 10 1 191-24-2 330 67 10000 Benzo(g,h,i)perylene

<sup>\*</sup>Specific quantitation limits are highly matrix dependent. The quantitation limits listed herein are provided for guidance and may not always be achievable.

<sup>\*\*</sup>Quantitation limits and J-values listed for soil/sediment are based on wet weight. The quantitation limits calculated by the laboratory for soil/sediment, calculated on dry-weight bases as required by the contract, will be higher.

<sup>&</sup>lt;sup>a</sup>The J-value for the medium-level soil analysis can be determined by multiplying the low-level soil J-value by a factor of 30.3.

J-values are evaluated annually and subject to change.

Section No. 9 Revision No. 5 Date: 05/12/97 Page 8 of 9

Target Compound List (TCL) and Contract Required quantitation Limits (CRQL)\* Low Soil (µg/kg) Water (µg/L) Quant. Quant. Pesticides/PCBs **CAS Number** Limit J-Value Limt\*\* J-Value\*\* alpha-BHC 319-84-6 0.05 0.0033 1.7 0.1 beta-BHC 319-85-7 0.05 0.0055 1.7 0.4 0.05 1.7 0.1 delta-BHC 0.013 319-86-8 gamma-BHC (Lindane) 58-89-9 0.05 0.0027 1.7 0.4 Heptachlor 76-44-8 0.05 0.002 1.7 0.1 1.7 Aldrin 309-00-2 0.05 0.002 0.1 Heptachlor epoxide 1024-57-3 0.05 0.0031 1.7 0.1 Endosulfan I 0.05 1.7 959-98-8 0.02 0.1 0.10 3.3 0.1 Dieldrin 60-57-1 0.015 4.4'-DDE 72-55-9 0.10 0.013 3.3 0.1 72-20-8 0.10 3.3 0.2 Endrin 0.018 0.10 3.3 Endosulfan II 33213-65-9 0.0084 0.3 4.4'-DDD 0.10 0.029 3.3 0.4 72-54-8 Endosulfan sulfate 1031-07-8 0.10 3.3 0.2 0.021 4,4'-DDT 50-29-3 0.10 0.015 3.3 0.1 17.0 Methoxychlor 72-43-5 0.50 0.12 1. Endrin ketone 53494-70-5 0.10 0.018 3.3 0.2 Endrin aldehyde 7421-36-3 0.10 0.018 3.3 0.2 5103-71-9 0.05 1.7 alpha-Chlordane 0.0028 0.1 gamma-Chlordane 5103-74-2 0.05 0.0031 1.7 0.1 170.0 Toxaphene 8001-35-2 5.0 0.2 10. Aroclor-1016 12674-11-2 1.0 0.14 33.0 3. Aroclor-1221 0.24 67.0 11104-28-2 2.0 8. Aroclor-1232 11141-16-5 1.0 0.2 33.0 5. Aroclor-1242 53469-21-9 1.0 0.51 33.0 6. Aroclor-1248 12672-29-6 1.0 33.0 0.16 3. Aroclor-1254 11097-69-1 1.0 0.04 33.0 4. Aroclor-1260 11096-82-5 1.0 0.15 33.0 2.

<sup>\*</sup>Specific quantitation limits are highly matrix dependent. The quantitation limits listed herein are provided for guidance and may not always be achievable.

<sup>\*\*</sup>Quantitation limits and J-values listed for soil/sediment are based on wet weight. The quantitation limits calculated by the laboratory for soil/sediment, calculated on dry-weight bases as required by the contract, will be higher.

J-values are evaluated annually and are subject to change.

Page 9 of 9

Inorganic Target /				
Analyte	Contract Required Detection Limit** (μg/L)			
Aluminum	200			
Antimony	60			
Arsenic	10*			
Barium	200			
Beryllium	5			
Cadmium	5			
Calcium	5000			
Chromium	10			
Cobalt	50			
Copper	25			
Iron	100			
Lead	3*			
Magnesium	5000			
Manganese	15			
Mercury	0.2			
Nickel	40			
Potassium	5000			
Selenium	5*			
Silver	10			
Sodium	5000			
Thallium	10*			
Vanadium	50			
Zinc	20			
Cyanide	10			

<sup>\*</sup>Graphite furnace or Trace ICP required.

Instrument Detection Limits (IDLs) are available upon request. IDLs are instrument specific and updated quarterly.

<sup>\*\*</sup>The CRDLs are the minimum levels of detection acceptable under the CLP SOW procedures. The detection limits for samples may be considerably higher depending on the sample matrix.

Section No. 10 Revision No. 2 Date: 07/25/95 Page 1 of 6

# 10. Data Reduction, Validation, and Reporting

Raw analytical data generated in the laboratories is collected on printouts from the instruments and associated data system or manually in bound notebooks. Analysts review data as it is generated to determine that the instruments are performing within specifications. This review includes calibration checks, surrogate recoveries, blank checks, retention time reproducibility, and other QC checks described in Sections No. 8 and No. 11. If any problems are noted during the analytical run, corrective action is taken and documented.

Each analytical run is reviewed by a chemist for completeness and accuracy prior to interpretation and data reduction. The following calculations are used to reduce raw data to reportable results.

GC/MS calculation used by the data system to determine concentration in extract for **semivolatiles** or in the sample itself for **volatiles**:

$$Q= (Ax) (Is)/(Als) (RRF) (Vi)$$

Where:

Ax = Peak area

Als = Internal standard peak area

Is = Amount of internal standard injected (ng)

RRF = Relative response factor

Vi = Volume of extract injected (L) or volume sample purged (mL)

Section No. 10 Revision No. 2 Date: 07/25/95 Page 2 of 6

The extract concentration is further reduced by considering the initial sample weight or volume and the final extract volume:

Concentration = 
$$(Q)$$
  $(D)$   $(F)$   $(1000)/(I)$ 

Where:

Q = Concentration determined by the data system (mg/L)

D = Dilution factor if needed

F = Final extract volume (mL)

I = Initial sample weight (grams) or volume (mL)

Results are reported in  $\mu$ g/L for water samples and  $\mu$ g/kg for solid samples. Soil samples are reported on an as received and on a dry-weight basis. The results are reported on Form I shown in Appendix A.

The results for the **pesticides/PCBs** analysis are calculated using the following equation:

Concentration = 
$$(Ax)$$
 (Is)  $(Vt)$   $(DF)/(As)$   $(Vi)$   $(Vs)$ 

Where:

Ax = Peak height for the parameter being measured

Is = Amount of standard injected (ng)

Vt = Volume of total extract ( $\mu$ L)

DF = Dilution factor, if needed

As = Peak height for the external standard

Vi = Volume of extract injected (µL)

Vs = Volume (mL) or weight (gm) of sample extracted

Section No. 10 Revision No. 2 Date: 07/25/95 Page 3 of 6

Results are reported as  $\mu g/L$  for water samples and mg/kg for solid samples. Soil samples are reported on a dry-weight basis. Results are reported on Form I shown in Appendix A.

The results for inorganic analyses are calculated using the following equation:

Concentration = (A)(D)(E)(1000)/(F)

Where:

A = The concentration determined by AA, ICP, or FTIR using calibration data programmed into the instrument (mg/L)

D = Dilution factor if needed

E = Final extract volume (mL)

F = Initial sample volume (mL) or weight (gm)

Results are usually reported in  $\mu g/L$  for water samples and in mg/kg for solid samples. Soil samples are reported on a dry-weight basis. The results are reported on Form I shown in Appendix A.

The principle criteria used to validate data will be the acceptance criteria described in Sections No. 8 and 11 and protocols specified in laboratory SOPs. Following review, interpretation, and data reduction by the analyst, data is transferred to the laboratory sample management system either by direct data upload from the analytical data system or manually. This system stores client information, sample results, and QC results. A security system is in place to control access of laboratory personnel and to provide an audit trail for information changes. The data is again reviewed by the group leader or another analyst whose function is to provide an independent review and verified on the sample management system. The person performing the verification step reviews all data including quality control information prior to verifying the data. Any errors identified and corrected during the review process are documented and addressed with appropriate personnel to ensure generation of quality data. If data package deliverables have been requested, the laboratory will complete the appropriate forms (see Appendix A) summarizing the quality control information, and transfer copies of all raw data (instrument printouts, spectra, chromatograms, laboratory notebooks, etc.) to the

Section No. 10 Revision No. 2 Date: 07/25/95 Page 4 of 6

Data Packages Group. This group will combine the information from the various analytical groups and the analytical reports from the laboratory sample management system into one package in the client requested format. This package is reviewed by the Quality Assurance Department for conformance with SOPs and to ensure that all QC goals have been met. Any analytical problems are discussed in the case narrative, which is also included with the data package deliverables.

The validation of the data by the Quality Assurance Department includes spot checking raw data versus the final report, checking that all pertinent raw data is included and does refer to the samples analyzed, review of all QC results for conformance with the method, and review of the case narrative for description of any unusual occurrences during analysis. This validation is performed using techniques similar to those used by the Sample Management Office for the USEPA's Contract Laboratory Program. The validation performed by the laboratory does not address usability of the data, which usually requires some knowledge of the site. The laboratory will make every attempt to meet the requirements of this QAPP, thus reducing the need to assess usability of the data.

The laboratory sample management system is programmed to accept and track the results of quality control samples including blanks, surrogates, recoveries, duplicates, controls, and reference materials. The computer is programmed with the acceptance criteria for each type of QC sample and will display an out-of-spec message if the data is not within specifications. All data outside of specifications appears on a report to the Quality Assurance Department on the next working day. These are reviewed by the Quality Assurance Department for severity of the problems and trends in the data. The reports are then sent to the analytical groups for the purpose of documenting the corrective action taken. The sample management system also produces control charts and has searching capabilities to aid in data review. The flow of data from the time the samples enter the laboratory until the data is reported are summarized in Table 10-1.

Any data recorded manually will be collected in bound notebooks. All entries will be in ink, with no erasures or white-out being permitted. Any changes in data will be made using a single line to avoid obliteration of the original entry and will be dated and signed. Any data resulting from instrument printouts will be dated and will contain the signature and/or identification of the analyst responsible for its

Section No. 10 Revision No. 2 Date: 07/25/95 Page 5 of 6

generation. After copies of the data are incorporated into the data package deliverables, the originals will be stored in locked archives at the laboratory for a period of 7 years.

Project files will be created per client/project and will contain chain-of-custody records, analysis requirements, and laboratory acknowledgments which document samples received, laboratory sample number assignment, and analysis requested. Raw data is filed per batch number assignment and laboratory sample number which correlates to the sample receipt documents. When the project is complete, all documentation is archived in a limited access area and retained for 5 years.

Section No. 10 Revision No. 2 Date: 07/25/95 Page 6 of 6

Table 10-1						
Sample and Data Routing at Lancaster Laboratories						
Action	Personnel Involved					
Sample received at Lancaster Labs	Sample Administration					
Sample is entered onto sample management system (lab ID number assigned, analyses scheduled, chain of custody started, storage location assigned)	Sample Administration					
Sample stored in assigned location (refrigerator, freezer, etc.)	Sample Support					
Acknowledgment sent to client	Sample Administration					
Removed from storage for analysis; necessary aliquot taken and sample returned to storage	Technical Personnel					
Analysis is performed according to selected analytical method; raw data recorded, reviewed, and transferred to computer by chemist or technician*	Technical Personnel					
Computer performs calculations as programmed according to methods	Data Processing					
Chemist or supervisor verifies raw data	Technical Personnel					
Data package deliverables are assembled	Data Package Group					
Data packages are reviewed prior to mailing	Quality Assurance Dept. Laboratory Management					

<sup>\*</sup>Analyses requiring the chemist's interpretation may involve manual data reduction prior to entry onto the computer.

Section No. 11 Revision No. 3 Date: 05/12/97 Page 1 of 9

#### 11. Internal Quality Control Checks

The particular types and frequencies of quality control checks analyzed with each sample are defined in USEPA CLP SOW OLM03.2 for organics and in CLP SOW ILM04.0 for inorganics or most recent revision, along with the limits of acceptance or rejection. The quality control checks routinely performed during sample analysis include surrogates, matrix spikes, duplicates, blanks, and internal standards. In addition to these checks, inorganic analyses employ serial dilutions, interference check samples, and laboratory control samples.

<u>Surrogates</u> (used for organic analysis only) - Each sample, matrix spike, matrix spike duplicate, and blank are spiked with surrogate compounds prior to purging and extraction in order to monitor preparation and analysis. Surrogates are used to evaluate analytical efficiency by measuring recovery.

<u>Matrix Spikes</u> - A matrix (soil or water) is spiked with known quantities of specific compounds and subjected to the entire analytical procedure in order to indicate the appropriateness of the method for the matrix by measuring recovery.

<u>Duplicates</u> (matrix spike duplicate - organics; duplicate - inorganics) - A second aliquot of a matrix/sample is analyzed at the same time as the original sample in order to determine the precision of the method. Recovery of the original compared to the duplicate is expressed as relative percent differences (RPD).

Blanks (method, storage, instrument) - Blanks are an analytical control consisting of a volume of deionized, distilled laboratory water for water samples and all storage and instrument blanks, or a purified solid matrix for soil/sediment samples. (Metals use a digested reagent blank with soils.) They are treated with the same reagents, internal standards, and surrogate standards and carried through the entire analytical procedure. The blank is used to define the level of laboratory background contamination.

Section No. 11 Revision No. 3 Date: 05/12/97 Page 2 of 9

Internal Standards (used for GC/MS analysis) - Internal standards are compounds added to every standard, blank matrix, spike, matrix spike duplicate, and sample at a known concentration, prior to analysis. Comparison of the peak areas of the internal standards are used for internal standard quantitation as well as to determine when changes in the instrument response will adversely affect quantification of target compounds.

<u>Serial Dilutions</u> (used for inorganic ICP analysis) - If the analyte concentration is sufficiently high (≥50× IDL) an analysis of a five-fold dilution must agree within 10% of the original determination. If the dilution analysis is not within 10%, a chemical or physical interference effect should be suspected.

Interference Check Sample (used for inorganic ICP analysis) - To verify interelement and background correction factors a solution containing both interfering and analyte elements of known concentration is analyzed at the beginning and end of each analysis run and per 20 samples.

<u>Laboratory Control Samples</u> (used for inorganic analysis) - Aqueous and solid control samples of known composition are analyzed using the same sample preparation, reagents, and analytical methods employed for the sample. LCS recovery must fall within established control limits.

The results of quality control samples are entered into the computer along with sample results. The computer is programmed to compare the individual values with the acceptance limits. If the results are not within the acceptance criteria, appropriate corrective action is taken where necessary. Management is kept informed by daily reports of QC outliers generated by the computerized system. Monthly reports on results of all QC analyses showing mean and standard deviation will indicate trends or method bias. Control charts are plotted via computer and may be accessed at any time by all analysts.

The tables that follow show the types and frequency of QC performed, along with the acceptance limits and corrective action.

Section No. 11 Revision No. 3 Date: 05/12/97 Page 3 of 9

#### Table 11-1

# Quality Control GC/MS Volatiles

		Limita/9/	7 <b>3</b>	T
Туре	Acceptance WATERS	SOILS	Frequency	Corrective Action
Surrogates:			Each sample, MS,	Reanalyze sample if outside
Taluana do	00 440	04 420	MSD, and blank	limits; if reanalysis confirms
Toluene-d8	88 - 110	84 - 138		original, document on report
Bromofluorobenzene	86 - 115	59 - 113		and/or case narrative
1,2-Dichloroethane-d4	76 - 114	70 - 121		ļ
Matrix Spikes:			Each group (≤20) of samples per	Advisory Only
1,1-Dichloroethene	61 - 145	59 - 172	matrix/level	Evaluated by analyst in
Trichloroethene	71 - 120	62 - 137		relationship to other QC
Benzene	76 - 127	66 - 142		results
Toluene	76 - 125	59 - 139		
Chlorobenzene	75 - 130	60 - 133		
Matrix Spike Duplicates		·	Each group (≤20) of	Advisory Only
(RPD):			samples per	
,			matrix/level	Evaluated by analyst in
1,1-Dichloroethene	14	22		relationship to other QC
Trichloroethene	14	24		results
Benzene	11	21		
Toluene	13	21		]
Chlorobenzene	13	21		j
Blanks:	<(2.5×) CRQL	for	Once for each	Reanalyze blank and
	methylene chlo		12-hour time period	associated samples if blank
1	induity to the		)	outside limits
•	<(5×) CRQL for acetone			
	and 2-butanon	1 , ,		
	and 2-butanon	C		
	<crql all="" for="" of<="" td=""><td>other</td><td></td><td></td></crql>	other		
	compounds	J		
Internal Standards:	-50% to +100%	6 of internal	Each sample, MS,	Reanalyze samples; if
mitomai Otamai US.	standard area		MSD, and blank	reanalysis confirms original,
Bromochloromethane	STD	01 12-110ui	MOD, and blank	document on report or case
1,4-Difluorobenzene	310			narrative
•	DT Change <2	0.000		Hallative
Chlorobenzene-d5	RT Change ≤3	U 360.	<u> </u>	

Some of the CLP recovery limits are advisory limits only. If in the opinion of the analyst, a problem exists other than matrix, the samples will be reanalyzed.

Section No. 11 Revision No. 3 Date: 05/12/97 Page 4 of 9

### Table 11-2

# Quality Control GC/MS Semivolatiles

GC/MS Semivolatiles						
Туре	Acceptance WATERS	Limits (%) SOILS	Frequency	Corrective Action		
Surrogate:			Each sample, MS,	Repeat analysis if more than		
			MSD, and blank	one surrogate out per fraction		
Nitrobenzene-d5	35 - 114	23 - 120		(acid/base) or any recovery		
2-Fluorobiphenyl	43 - 116	30 - 115		<10%; if reanalysis confirms		
Terphenyl-d14	33 - 141	18 - 137		originals, document on report		
Phenol-d5	10 - 110	24 - 113		and/or case narrative		
2-Fluorophenol	21 - 110	25 - 121		1		
2,4,6-Tribromophenol	10 - 123	19 - 122				
2-Chlorophenol-d4	33 - 110	20 - 130	(Advisory)	•.,		
1,2-dichlorobenzene-d4	16 - 110	20 - 130	(Advisory)			
Matrix Spikes:			Each group (≤20) of	Advisory Only		
•			samples per			
Phenol	12 - 110	26 - 90	matrix/level	Evaluated by analyst in		
2-Chlorophenol	27 - 123	25 - 102		relationship to other QC		
1,4-Dichlorobenzene	36 - 97	28 - 104		results		
N-Nitroso-di-n-propylamine	41 - 116	41 - 126		•		
1,2,4-Trichlorobenzene	39 - 98	38 - 107		1		
4-Chloro-3-methylphenol	23 - 97	26 - 103				
Acenaphthene	46 - 118	31 - 137		\ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \		
4-Nitrophenol	10 - 80	11 - 114				
2,4-Dinitrotoluene	24 - 96	28 - 89				
Pentachlorophenol	9 - 103	17 - 109		1		
Pyrene	26 - 127	35 - 142				
Matrix Spike Duplicates (RPD):	-		Each group (≤20) of	Advisory Only		
,	ţ		samples per			
Phenol	42	35	matrix/level	Evaluated by analyst in		
2-Chlorophenol	40	50		relationship to other QC		
1,4-Dichlorobenzene	28	27		results		
N-Nitroso-di-n-propylamine	38	38		<u>'</u>		
1,2,4-Trichlorobenzene	28	23				
4-Chloro-3-methylphenol	42	33				
Acenaphthene	31	19				
4-Nitrophenol	50	50				
2,4-Dinitrotoluene	38	47		į		
Pentachlorophenol	50	47				
Pyrene	31	36				
Blanks:	<(5×) CRQL fo		Once per case or	Reextract and reanalyze blank		
	phthalate este		group (≤20) of	and associated samples		
			samples, each	1,		
	<crql all<="" for="" td=""><td>other TCL</td><td>matrix, level,</td><td></td></crql>	other TCL	matrix, level,			
	compounds	: - <del></del>	instrument	1		
	33504.140		1.13d dillone	<u> </u>		

Section No. 11 Revision No. 3 Date: 05/12/97 Page 5 of 9

**Table 11-2** 

Quality Control GC/MS Semivolatiles

	Acceptance Limits	(%)	
Туре	WATERS SO	ILS Frequency	Corrective Action
Internal Standards:	-50% to +100% of inte standard area of 12-ho	1 ' '	Reanalyze samples; if reanalysis confirms original,
1,4-Dichlorobenzene-d4 Naphthalene-d8	STD		document on report and/or case narrative
Acenaphthene-d10 Phenanthrene-d10	RT change ≤30 sec.		
Chrysene-d12 Perylene-d12			

Some of the CLP recovery limits are advisory limits only. If in the opinion of the analyst, a problem exists other than matrix, the sample will be reanalyzed.

Section No. 11 Revision No. 3 Date: 05/12/97 Page 6 of 9

#### **Table 11-3**

### Quality Control Pesticides/PCBs

		Pesticides/PCI	Bs	
		e Limits (%)		
Туре	WATERS	SOILS	Frequency	Corrective Action
Surrogates:			Added to each sample, MS/MSD,	Advisory Only for Samples
Tetrachloro-m-xylene	30 - 150	30 - 150	and blank during the	For Blank, reinject; if still out
Decachlorobiphenyl	30 - 150	30 - 150	extraction phase	reextract and reanalyze blank and associated samples
Matrix Spikes:			Each extraction group (≤20) of	Advisory Only
gamma-BHC (Lindane)	56 - 123	46 - 127	samples per	Evaluated by analyst in
Heptachlor	40 - 131	35 - 130	matrix/level	relationship to other QC
Aldrin	40 - 120	34 - 132		results
Dieldrin	52 - 126	31 - 134		
Endrin	56 - 121	42 - 139		
4,4'-DDT	38 - 127	23 - 134		
Matrix Spike Duplicates (RPD):			Each group (≤20) of samples per	Advisory Only
gamma-BHC (Lindane)	15	50	matrix/level	Evaluated by analyst in
Heptachlor	20	31		relationship to other QC
Aldrin	22	43		results
Dieldrin	18	38		
Endrin	21	45		
4,4'-DDT	27	50		
Blanks:	<pre><crql an<="" for="" pre=""></crql></pre>		Once per case or	Inject a hexane or solvent
	compounds	•	extraction group	blank first to be sure the
	, i		(≤20) of samples,	analytical system is clean then
	<(.5×) the CR	QL for	each matrix, level,	reinject the blank itself. If the
	instrument bla		instrument	reinjected blank is acceptable,
				any samples extracted with
				this blank should be reinjected
				if they, too, contain the analyte
				which was contaminating the
				blank. If the reinjected blank
				is unacceptable, any affected
				samples must be reprepped.

Some of the CLP recovery limits are advisory limits only. If in the opinion of the analyst, a problem exists other than matrix, the sample will be reanalyzed.

Section No. 11 Revision No. 3 Date: 05/12/97 Page 7 of 9

### Table 11-4

# Quality Control Inorganics

	Inorganics		<b>-</b>
Туре	Acceptance Limits (%) WATERS SOILS	Frequency	Corrective Action
Matrix Spikes:	75% to 125% except where sample conc. exceeds spike conc. by ≥4×	Each group of samples of similar matrix/level (≤20) each method	Analyze post-digestion spike sample
Duplicates (RPD):	±20% RPD for sample values ≥5× CRDL	Each group of samples of similar matrix/level (≤20) each method	Flag the data
Blanks: Initial Calibration (ICB) Continuing Calibration (CCB)	≤CRQL	Each wavelength immediately after calibration verification at 10% frequency or every 2 hours (beginning and end of run min.)	Correct problem, recalibrate, and rerun
Preparation Blank	≤CRDL >CRDL then lowest conc. in sample must be 10× blk. conc. or <crdl< p=""></crdl<>	Each SDG or batch (≤20 samples)	Redigest and reanalyze blank and associated samples if sample result <10× blank result
Serial Dilutions:	Within ±10% of the original determination	Each group of (≤ 20) of similar matrix/level	Flag the data
Interference Check Sample:	Solution A - ± (2×) CRDL of the true value for analytes with CRDLs of ≤10 μg/L Solution AB - ±20% of the true value for the analytes	Each wavelength after Initial Calibration Verification at beginning and end of the run and per 20 samples	Recalibrate the instrument
Laboratory Control Sample:	Aqueous 80% to 120% (except Ag and Sb)  Solids see Table 11-15	Each SDG or batch (≤20 samples), each method	Redigest and reanalyze LCS and associated samples
Post Digestion Spike:	85% to 115%	When matrix spikes are outside 75% to 125% range (not performed on GFAA analyses)	Flag the data
Analytical Spike:	85% to 115%	Every GFAA determination	See Figure 11-1

Section No. 11 Revision No. 3 Date: 05/12/97

Page 8 of 9

**Table 11-5** 



### Certification

### PriorityPollutnT™/CLP Inorganic Soils

Quality Control Standards

Catalog Nº PPS-46

Lot Nº 229

Parameter	Certified Value	Performance Acceptance Limits™
TRACE METALS PriorityPollutn™	mg/Kg	mg/Kg
(Catalog No 540)		
aluminum	4590	2280 - 7590
antimony	39.8	8.37 - 119
arsenic	75.4	37.1 - 112
barlum	106	74.3 - 139
beryllium	51.0	11.7 - 90.3
boron	94.1	26.9 - 161
cadmium	45.4	11.9 - 79.0
calcium	1290	875 - 1750
chromium	71.0	38.0 - 100
cobalt	49.5	29.8 - 70.5
copper	112	63.9 - 162
iron	9160	5560 <b>- 13000</b>
lead	53.5	28.1 - 75.9
magnesium	1160	691 - 1670
manganese	154	107 - 208
mercury	1.50	0.389 - 2.35
molybdenum	47.4	29.2 - 70.2
nickel	39.4	21.5 - 57.5
potasskim	1420	880 - 1870
selenium	72.3	37.8 - 108
siver	116	53.2 - 170
sodium	198	111 - 287
strontium	109	46.3 - 173
thallum	40.0	20.0 - 60.0
tia	102	35.9 - 168
titanium	230	60.0 - 400
vanadium	€5.9	32.0 - 88.9
zinc	134	72.2 - 199
CYANIDE PriorityPollutnT™	mg/Kg	mg/Kg
(Catalog No 541)	- •	
total cyanide	323	123 - 559

The Trace Metals Certified Values are equal to the mean recoveries for each parameter as determined in an interlaboratory round robin study. The standard was digested using Method 3050, SW-846 and the digest analyzed by ICP and stomic absorption spectroscopy.

The Cyanide Certified Value is equal to the mean recovery as determined in an interlaboratory round robin study. The standard was distilled and analyzed following the procedure outlined in Method 9010, SW-846.

The Performance Acceptance Limits (PALs™) are listed as guidelines for acceptable analytical results given the limitations of the USEPA methodologies commonly used to determine these parameters and closely approximate the 95% confidence interval. The PALs™ are based on data generated by your peer laboratories in ERA's InterLaB™ program using the same samples you are analyzing and data from USEPA methods, WP, WS and CLP interlaboratory studies. If your result falls outside of the PALs™, ERA recommends that you investigate potential sources of error in your preparation and/or analytical procedures. For further technical assistance, call ERA at 1-800-372-0122.


For users of internal standards, ERA has determined that scandium is present in this soil at 1.66 mg/Kg and that yttrium is present at 9.43 mg/Kg.

\*Each lot of standards will have different certified values and the advisory range will be adjusted accordingly.

Section No. 11 Revision No. 3 Date: 05/12/97 Page 9 of 9

Figure 11-1

FIGURE 1. FURNACE ATOMIC ABSORPTION ANALYSIS SCHEME



Section No. 12 Revision No. 3 Date: 10/17/96 Page 1 of 18

#### 12. <u>Performance and System Audits</u>

System audits are conducted on each department at Lancaster Laboratories by members of the Quality Assurance Department. The audits include checks on methodology, reagent preparation, equipment calibration and maintenance, quality control results, and training of personnel. The results of the audits and corrective action, where necessary, are communicated to laboratory personnel and management by means of a written report. Audits by outside organizations including clients, regulatory personnel, and the USEPA are permitted by arrangement with the Quality Assurance Department.

The Quality Assurance Department reviews summaries of the quality control data entered onto the computerized sample management system by analysts. Control charts and statistics are reviewed for trends which may indicate problems with the analytical data. In this way, small problems are identified before they have any significant impact on laboratory results.

Performance audits consist of both intralaboratory and interlaboratory check samples. QC samples from commercial suppliers are analyzed quarterly to assess laboratory accuracy including a double blind program. The Laboratory also participates in a number of interlaboratory performance evaluation studies which involve analysis of samples with concentrations of analytes that are known to the sponsoring organization, but unknown to the laboratory. Inorganics, pesticide/herbicides, trihalomethanes, volatile organic compounds, semivolatile organic compounds, and traditional wet chemistry analyses are analyzed by Lancaster Labs for studies conducted by the USEPA and the New York Department of Health. Lancaster Labs has participated in the USEPA Contract Laboratory Program which provides laboratory analysis in support of the Superfund program. Part of maintaining this contract includes analysis of quarterly blind samples. Representative results from some of these studies are attached to this section.

Section No. 12 Revision No. 3 Date: 10/17/96 Page 2 of 18

# LANCASTER LABORATORIES Account# 7174662301 LOW LANC LANCASTER PA GCL

Performance Evaluation Report USEPA Water Supply Study WS037

Page: FF005 Page: 1 Date: 065EP96

Participant	ID: PA00009 Type: OTHER			Requesting (	Office: U1
	Sample Number	Reported Value	True Value⇒	Acceptance limits	
TRACE ME	TAIS TN M	ICROGRAMS PI	PR TTTFR:		
01-ARSENIC					
002-BARIUM	001	049.0	49.3	41.9- 56.3	Accept.
	002	0771.	773	657- 889	Accept.
03-CADHIUM	001	C10.2	10.2	8.16- 12.2	Accept.
04-CHROMIUM					•
05-LEAD	001	071.5	72.9	62- 83.8	Accept.
AC MEDANA	001	013.2	13.8	9.66- 17.9	Accept.
06-HERCURY	001	07.70	e.16	5.71- 10.6	Accept.
07-SELENIUM	001	051.3	57.9	46.3- 69.5	Accest.
91-COPPER		037.03		40.3- 07.3	ACCE L.
40-ANTIMONY	CC1	054.0	55.7	50.1- 61.3	Accept.
- WILLIAM	002	021.4	18.0	12.6- 23.4	Accept.
.41-BERTLLIU	M 001	03.27	4.26	3.62- 4.9	Not Accep
42-NICKEL					
43-THALLIUM	001	055.9	55.0	46.8- 63.3	Accept.
	002	C2.40	2.38	1.67- 3.69	Accept.
26-BORON	002	0953.	929	876- 1030	Acc€ŗt.
36-MANGANESI	2				·
37-noly bde ni	C01	047.8	48.1	43- 51.4	Accest.
	002	053.1	54.0	42.6- 65.4	Accept.
39-ZINC	001	0588.	600	536+ 652	Accept.
NITBATE/!   NITBATE = 0		LUCRIDE IN 1	HILLIGRAMS .	PPR LITER:	
	001	08.45	8.30	7.47- 9.13	Accept.
92-WITEITE	AS N 001	0.493	0.502	0.427~0.577	Accept.
61-ORTHOPHO	STHATE AS	P			•
	001	C1.11	1.10	0.957- 1.21	Accept.
	IDES IN M	ICROGRAMS PI	ER LITER:		
11-ENDEIN	001	0.301	0.231	0.162- 0.3	Not Accep
	001	04341	0.431	90104 903	HOL FCC

Section No. 12 Revision No. 3 Date: 10/17/96 Page 3 of 18

Performance Evaluation Report USEPA Water Supply Study WS037 Report: PEOOS Page: ? Date: O65FP96

Participant ID: PA00009		)9 T	pe: OTHER	Requesting Office: UT	
	Sample Number	Reported Value	True Value‡	Acceptance Limits	Performance Evaluation
012-LINDANE					
013-HETROXYCE		C.353	0.381	0.21-C.552	Accept.
014-TOXAPHENE	. 001	015.5	18.5	10.2- 26.8	Accept.
093-ALACHION	002	07.04	8.81	4.85- 12.8	Accept.
O) J A LACILLON	005	05.36	4.87	2.68- 7.06	Accept.
094-ATRAZINE	005	07.30	6.80	3.74- 9.86	Accept.
095-HEPTACHLO		A 267	0 543	0 31-0 816	Accort
096-HEPTACHLO	004 F POXIDE	C.367	0.563	C.31-0.816	yccutt.
	004	C.4C6	0.403	0.222-0.584	Accept.
097-CHLCHDAHE	(TOTAL) 003	02.39	4.44	2.44- 6.44	Not Accept.
113-SIMAZINE	005	<b>C6.30</b>	5.56	1.04- 9.77	Accest.
172-HEXACHLOR		00.30	3.30	1.04- 3.77	nece ; c ·
844	004	0.618	0.806	0.323- 1.14	Accept.
241-METOLACHL 242-METRIBUZI	006	021.7	19.4	7.87- 29.5	Accept.
	006	C14.9	14.1	D.L 22.4	Accept.
243-PROMETON	006	023.0	18.8	6.48- 28.3	Accept.
256-ALDRIN	004	0.433	0.567	0.186-0.725	Accept.
257-BUTICHLOR	006	C22.3	20.5	5.93- 31.3	Accept.
258-DIELDRIN					•
259-PROPACHLO	004 H	0.554	0.530	0.358-0.708	Accept.
	004	01.16	1.20	C.566- 1.86	Accept.
CABBANATE 098-ALDICARB	S IN MICH	OGRAMS PER	LITER:		
U) REDICKIO	C01	C36.4	34.3	24.3- 41.4	Accept.
099-ALDICARB	SULFONE CO1	034.1	32.1	28.7- 40.1	Accept.
100-ALCICARB			36.	20.7- 40.1	wcct.
	001	027.6	25.9	20.3- 33	Accept.
101-CAREOPURA		040 "	4.0 0	26.0.26.2	
114-OXARYL (V	001 V(AT4)	042.4	48.9	26.9- 76.9	Accept.
TTA CHARITY (A	001	044.0	46.4	36.3- 54.9	Accept.

Section No. 12 Revision No. 3 Date: 10/17/96 Page 4 of 18

Performance Evaluation Report USEPA Water Supply Study WS037

Report: FE005 Page: Date: 065FP96

Participant 1	IC: PAOOO	09 τη	pe: OTHER	Fequesting C	ffice: UT
	Sample Number	Reported Value		Acceptance Limits	Performance Evaluation
245-METHONYL					
	001	060.0	60.7	49.4- 68.4	Accept.
	S IN MIC	ROGRAMS PER	LITER:		
015-2,4-0	001	013.0	14.9	7.45- 22.4	Accept.
016-2,4,5-TP			2 / 2 /		•
	001	09.53	11.0	5.9- 17.7	Accept.
102-PENTACHLO		05.07	6.59	3.3- 9.89	Accept.
115-DALAPON	001	03.07	0.39	7. 7- 2.65	ACCO; C.
	002	047.1	56.4	D.I 94.8	Accept.
116-DINOSEB					
117 DYCLODLA	002	014.2	18.6	C.652- 29.6	yccelt.
117-PICLORAN	002	017.9	23.3	D.L 34.8	Accept.
247-DICAMBA	002	017.7	23.3	5466	"CCC ; "
	002	031.5	30.4	2.98- 58.7	Accept.
POLYCHLOR 118-DECACHLOR			MICHOGRAMS		
	001	C.305	0.527	D.I 1.05	Accept.
		MS FER LITEE	<b>:</b>		
122-BENZO (A) P	CO1	C.754	C.937	0.115- 1.31	Accept.
		S IN MICROGE	RAMS PER LI	TER:	
134-DI (2-ETH			24.2	11 4 60 3	
136-DI (2-ETH	CO1	026.7	34.3	11.4- 52.3	Accept.
130-21 (1 111	001	016.6	21.3	6.98- 34.5	Accept.
MISCELLAN 137-Diquat	EOUS SOC	s in micro	GRAMS PER L	ITER:	
	C01	C3.43	8.41	2.05- 22.4	Accept.
138-endothall	001	C98.6	179	12- 312	Accest.
139-GLYFHOSAT		4,500	2.,,	11 311	
	001	0729.	780	630- 903	Accept.
TRINALOME 017-CHLOROFOR		N MICROGRAMS	FER LITER	•	
	001	024.1	22.3	17.8- 26.8	Accept.
018-B RONO PO RM	001	C18.9	18.6	14.9- 22.3	Accept.
					•

Section No. 12 Revision No. 3 Date: 10/17/96 Page 5 of 18

Performance Evaluation Report USEFA Water Supply Study WS037

Peport: PPOOS Page: . u Date: O65EP96

Participant	ED: PA000	09 T	ype: OTHER	Requesting (	office: UT
		Reported Value		Acceptance Limits	Performance Fvaluation
019-BROMODIC	HLOROMETH	ANE			
	001	C12.2	12.7	10.2- 15.2	Accept.
0 20 - CHL CROD	BRCHCHETH				
	001	015.3	14.2	11.4- 17	Accept.
021-TOTAL TI	HTSHOLAHIS				
	001	070.5	67.8	54.2- 81.4	Accept.
VOLATILI	ORGANIC	COMPCUNES I	N MICROGRAM	S PER LITFF:	
032-VINTL CH					
	001	015.5	14.8	8.88- 20.7	Accept.
034-1,1-DIC	LORGETRYL	ENE			
	001	018.3	16.5	13.2- 19.8	Accept.
035-1, 2-DICH	LCROETHAN				
	001	015.9	13.2	10.6- 15.8	Rot Accept.
036-1,1,1-7	ICHLOROET	H A N E	`_		
	COI	C11.9	10.3	8.24- 12.4	Accept.
037-CARBON T	ETRACHLOR:	IDE			
	001	014.5	12.7	10.2- 15.2	ACCERT.
038-TRICHLOR	OETHYLENE			•	•
	001	08.20	8.70	5.22- 12.2	Accept.
039-BENZENE					
	001	013.0	12.5	10- 15	.ccept.
040-TETBACHI	ORCETHTLE	NE		:	
	002	010.1	9.60	5.76- 13.4	Accept.
041-1,4-DICH	LOROBENZE	NE			
	001	C6.65	7.31	4.39- 1C.2	Accept.
042-T 1,2 DI	CHLOROETH	YLENE			
	002	C15.0	14.8	11.8- 17.8	Accept.
043-C 1,2 DI	CHLOROETH	TLENE			
	002	011.4	9.72	5.83- 13.6	Accept.
044-1,2 DICH	LO ROPEOPA	N P			
	002	015.4	14.2	11.4- 17	Accept.
045-1,2DIBRO	MC3CHLORO	PROFANE			
	004	0.274	0.286	0.172- 0.4	Accept.
046-ethtlene	DIBRONID	E (ECH)			
	004	C.151	C.138	0.0828-0.193	Accept.
047-TOLUENE					
	002	C5.74	5.70	3.42- 7.98	Accept.
048-RTHYLBEN	ZENE				-
	<b>002</b>	09.40	9.19	5.51- 12.9	Accept.
049-CHLOROBE	NZENE				
	002	08.42	8.31	4.99- 11.6	Accept.
053-STYPENE					
	002	07.60	7.40	4.44- 10.4	Accept.
054-1,2 DICE	LOROBENZE				•
-	002	014.3	14.5	11.6- 17.4	Accept.

Section No. 12 Revision No. 3 Date: 10/17/96 Page 6 of 18

Performance Evaluation Report USEPA Water Supply Study WS037

Report: FEOOS Page: 5 Date: CESEP96

Participant ID	: PA0000	9 т	ype: OTHER	Requesting	office: UT
	ample umber	Reported Value	True Valuo≑	Acceptance Limits	Ferformance Evaluation
055-DICHLOROME	THANE				
	001	09.48	8.41	5.05- 11.8	Accept.
056-1,1-DICHLO					
	003	C14 .8	13.6	11.2- 16.4	Accept.
061-1,1,2-TRIC					
	001	011.8	10.7	8.56- 12.8	Accept.
063-1,1,1,2TET					
	003	C17.0	15.3	12.8- 10	yccelt.
064-1,2,3-TRIC			4.20		
004 0 0 0 000	003	C8.32	8.29	1.53- 11	Accept.
076-1, 2, 4-TRIC			14. 3	11 " 17 3	
A33 4 3 3 FATC	002	014.7	14.3	11.4- 17.2	Accept.
077-1,2,3-TRIC			16 7	10.6- 21.2	Accept.
001 - 4 02 4 7 4 7 00 0	003	C15.7	16.7	10.4- 21.2	яссере.
081-HEXACHLORO			9.50	4.19- 14.4	Accept.
217 TARABAT VYIE	003	011.1	9.30	4.13- 14.4	Accete.
090-TOTAL XYLE	NES 002	015.3	12.9	10.3- 15.5	Accept.
152-C 1,3 DICH			12.7	10.3- 15	иссерс.
132-C 1,3 DICE	003	010.9	12.3	8.22- 14.3	Accept.
153-T 1,3 DICH			14.5	0.22- 14.5	жесере.
133-1 1,3 5101	003	016.4	17.5	11- 20.5	Accept.
			UCTS IN HIC	ROGRAMS PER LIT	FR:
157-DIBROMOACE					
	001	0.918	8.50	D.L 13.8	Acceşt.
158-DICHLOROAC	ETIC ACI				
	001	02.82	22.7	6.83- 30.3	Not Accept.
160-MCNCBROMOA					
	001	02.74	14.4	1.26- 21.4	Accept.
161-MONOCHLORO					
	001	C1.80	12.8	3.43- 21	Not Accept.
162-TRICHLOROA					
254	001	C3.77	32.3	5.47- 47.9	Not Accept.
250-BRONOCHLOR					
	001	02.36	19.8	3.19- 30.8	Rot Iccept.
INORGANIC	DISINFE	TION BY-PE	ODUCTS IN I	TICROGRAPS PER L	TTER:
193-BROMATE					
	002	05.68	4.56	D.L 29	Accept.
194-CHLORATE					
	001	<b>C92.1</b>	82.1	62.1- 100	Accert.
195-CHLCRITE					
	001	0165.	140	86.6- 213	Accept.
260-BRONIDE					
	002	0157.	140	113- 169	Accept.

Section No. 12 Revision No. 3 Date: 10/17/96 Page 7 of 18


Performance Evaluation Report USEPA Water Supply Study WS037

Report: 5E005 Page: 6 Date: CESED96

Participan	t ID: PAOOO	09 T	pe: OTHER	Requesting	Office: UT
	Sample Number	Reported Value		Acceptance Limits	Ferforwance Evaluation
HISCELI	LANECUS ANA	LYTES:			
		ORINE (MILLIC	GRAMS PER L	ITER)	
	001	02.64	2.20	2.03- 3.07	Accest.
023-TURBIDI					•
	001	61.40	1.54	1.26- 1.98	Accept.
024-TOTAL		RESIDUE (MILI			
	001	C254.		188- 434	Accept.
025-CALCIUM	· · · <del>-</del>	MG. CACU3/L)			
	001	0147.	144	137- 158	Accept.
026-PH-UNIT	rs .				·
	001	08.94	9.13	8.88- 9.31	Accept.
027-ALKALIN	HITY (MG. CA	CO3/L)			•
	001	028.6	27.4	25.7- 31.5	Accept.
029-SCDIUM	(HILLIG RAMS	PER LITER)			·
`	001	012.9	12.6	11.4- 13.7	Accest.
145-SULFATE	(MILLIGRAM	S PER LITER)	)		·
	001	C263.	280	253- 316	Accept.
146-TOTAL C	CYANIDE (MIL	LIGRAMS PER	LITER)		·
	0C1	0.337	0.380	0.285-0.475	Accept.
263-TOC					
	001	03.46	2.80	2.49- 3.24	Not Accept.

<sup>\*</sup> Based on gravimetric calculations, or a reference value when necessary.

Section No. 12 Revision No. 3 Date: 10/17/96 Page 8 of 18



# Performance Evaluation Report USEPA Water Pollution Study MP035

rage: 1 Date: 16APB96

Participant ID: PACCCC9			17 pe: 01 HER	naquesting Office: UT			
Sample Numter	Reported Value		Acceptance Limits	Varning Limits	rerformance Fvaluation		
TRACE METAL	S IN MICRO	I I Nemara	J , A				
001-ALUH1KUP	3.4.0	221	2/1 2/2	. 200			
C 1	310	321	261 - 362	?76- 367 1330- 1640	Accept.		
(2	1170	1500	1270- 1700	( ) yo = 16 aa	*Conbit		
UU2-ARSENIC Ul	192	101	167- 231	175- 223	Accept.		
(7	2 (	571	992- 676	£15- 6£3	*ccert.		
003-82RYLLIU4	. * 7	171	076	. 1 - ( :	"CUPIU"		
ן) נו	176	1.90	165- 209	170- 204	Accept.		
·. 3 C 2	526	541	nac- 597	495- EU3	Agrept.		
004-CADHIUK	_ Z (;	1	14 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1		* 14 T * 15 T *		
C1	50.4	52.6	nn.5- 60.7	46.5- 58.7	Accept.		
02	388	401	305- 959	359- 440	Accept.		
005-CCBAL1	31117		777 W 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	1,1 1	· · · · · · · · · · · · · · · · · · ·		
01	27.2	28.1	22.9- 32.6	20.1- 31.4	*ccent.		
Č 2	6(3	624	557- 686	574- 670	Accept.		
DO6-CHRONIUM	0.3		, , , , , , , , , , , , , , , , , , , ,	, .,,			
01	15.8	17.0	13- 20.5	13.9- 19.5	*ccept.		
C 2	8 5 C	88 C	767- 995	754- 958	iccert.		
OO7-COPPER		•					
<b>C1</b>	e3.2	06.7	75.5- 96.9	76.2- 94.2	Modert.		
C 2	358	170	334- 409	344- 399	Accept.		
NO8-180N					•		
01	41.0	30.4	18.8- 42.6	21.8- 39.6	Ct. for Frr.		
02	457	0.60	441- 519	1151- 509	≯ccept.		
DO9-HERCUFY			_		•		
01	3.36	3.10	2.03- 4.07	2.29- 3.01	Accept.		
(2	12.8	11.6	8.65- 14.7	9.41- 13.5	Accept.		
110-EANGANESE				1	•		
01	385	401	369- 441	178- 412	Accept.		
(2	865	P P 1	933~ 968	PSC- 951	Accept.		
D11-NICKEL				•	·		
G1	461	496	453- 560	466- 547	Accept.		
C 2	6 C 1	611	557- 6911	574- 680	Accept.		
012-LEAD					•		
<b>C1</b>	202	297	259- 334	269- 125	Accept.		
0.2	385	399	356- 446	767- 435	Accept.		
013-SELENIUF							
01	487	522	402- 615	029- 588	Accept.		
(2	898	ባንቦ	750- 1156	PC4- 1100	Accept.		
DIU- VANADIUP					•		
61	26.2	311	186- 239	192- 729	*ccept.		
(2	786	011	724- 886	745- 867	Accept.		
U15-21NC	, -	·			F		
01	71.9	71.9	62.7- 84.9	€ . c - N2.2	*ccept.		

Section No. 12 Revision No. 3 Date: 10/17/96 Page 9 of 18

Ferformance Evaluation Report Fage: 2
USEPA Water Pollution Study WPO35 Date: 16APR06

Nuster	Reported Value	True	Accent		Uar	. Laa	Dorfornaco
			Lie				Performance Fvaluation
D16-ANTIMENY							******
03	345	370	240-	45C	266-	473	Accept.
04	551	570	369-	692	1130-	651	Accept.
117-SILVER							
03	176	180	153-		160-	2 C C	Accept.
CH	337	340	298-	391	310-	JèU	/cc∈pt.
SULLIANT-81							
03	80.0	43.3	63.4-	99.1	67.9-	94 . €	Accept.
CH	356	165	30 <b>1</b> -	425	317-	410	Accept.
74-HOLTBDEAUK							
0 3	126	1 30	106-	151	112-	146	Accept.
04	309	310	257-		270-	345	Accent.
75-STRCNTIUM							•
03	3.5	3.55	2.56-	4.49	2.91-	4.23	teept.
04	94.0	96.0	70.8-		03.4-		Accept.
176-TXT # K X U F	. • •	- · ·	• •				· · ·
03	110	115	96.8-	130	191-	126	Accept.
(4	272	270	230-		239-		Accept.
019-PH-UNITS 63 04	4.34 5.57	5.50	5.46-	4.4 5.62	4.25- 5.48-		Accept.
320-SPEC. CCND. (	UFRCS/CH	AT 25 C)					
01	907.	916		983	6 it ∂ ~	9611	/ccert.
02	501.	586	536-	627	547-	616	Accept.
021-105 AT 180 C	-						
01	509.	553	326-	762	380-	700	tecept.
C 2	314.	311	226-	398	24P-	277	tecept.
22-TOTAL HARDNE	SS (AS CAC	:03)			1		
01	309.	3 3 0	302-		309-	351	Accept.
(2	57.2	101	90.8-	110	93.2-	108	Accept.
23-CALCIUM							
01	105	104	92.6-				*ccept.
02	6.63	6.39			5.79-		Accept.
24-MAGNES TUN							•
01	16.6	17.C	15.2-	17.3	15.7-	10.7	Accept.
02	20.7	20.6	18- 23		10.7-		Accent.
25-500101		<del>-</del>		-			-
01	14.8	10.2	13.1-	16.2	13.5-	15.8	Accept.
0 2	52.5	54.3	49.3-		50.5-		*ccept.
026-POTASSIUM	3	• •			. •		, J
01	21.4	21.0	111.15-	23.7	19.4-	23.1	Accept.
C 2	39.7	38.3	13.3-		34.3-		*ccent.
DZ7-TOTAL ALKALI					• • •	. , • •	~ · · · · ·
	-						
(1	21.0	20.0	17.4-	7 6 1	19.3-	211 1	Accept.

Section No. 12 Revision No. 3 Date: 10/17/96 Page 10 of 18

# Ferformance Evaluation Report USEPA Water Pollution Study MPC35

Tage: 3 Nate: 1615096

Farticipant ID:	PA00009		Type: CTHER	Requesti	ng Office: UT
			icceptance limits		Performance Evaluation
028-CHLORIDE					
<b>C1</b>	236.	241	224- 259	228- 254	Accept.
0 2	68.1	72.7	65.1- 79.8	67- 77.9	Accept.
C29-FIUORICE					
C 1	3.77	3.50	3.09- 3.8	3.10- 3.71	Ct. for Frr.
02	1.39	1.35	1.16- 1.53	1.21- 1.48	Accent.
030-SULFA1E					
01		18.9			
C 2	83.3	n6.4	72- 97	75.1- 93.9	Accept.
NUTRIERTS I NUTRIERTS I		MS/LITEF	i		
C1		19.0	15- 22.3	15.9- 21.5	Accept.
02	1.62		1.15- 2.09	1.26- 1.97	
C32-NITRATE-NIT	ROGEN				
<b>C1</b>	8.25	0.31	6.76- 9.69	7.11- 9.34	Accept.
02	.340	0.390	0.28-0.495	0.305-0.469	Accept.
033-ORTHOFHCSPH					·
C1	.054	.0560	0.0333-0.076	0.0384-0.071	Accept.
0 2	2.88	2.80	2.43- 3.19	2.52- 3.1	Accept.
034-KJELDAHL-HI	TREGEN				
03	.632		0.115- 1.12		
C 4	9.30	7.80	5.73- 9.64	6.2- 9.17	Cr. for Err.
035-TOTAL PHOSP		_			
(3		0.574	0.47-0.705	0.498-0.677	iccent.
C4	€.16	6.08	5.16- 7.2	5.41- 6.56	Accept.
DEMANDS IN 036-COD	MILLIGRAMS	ZETTER			
01	213.	236	109- 259	1981 250	Accept.
C 2	85.4	101	71.7- 120	77.3- 114	Accept.
037-TOC					•
C 1	91.2		78.5- 108	02.4- 105	Modert.
0 2	39.2	40.1	31.6- 47.4	?3.6- US.7	Accent.
C38-5-DAY 80D					
(1	152.		64.1- 218	166 - e i i	Accept.
02	58.7	62.5	29.5- 95.5	37.7- 87.7	tocent.
102-CARBONAC EOU					
C 1	167.	117	34.3- 199	55.6- 170	tampt.
C 2	65.4	51.6	20- 83.2	20.7- 75	tcaent.
FC8*5 1F M1		11 F F			
U42-PCR-AROCLOR				_	
<b>C1</b>	2.74	2.76	0.769- 4.3	1.17- 3.04	*ccert.
044-PCB-AROCLOR				•	
C 2	4.34	4.26	1.77- 6.04	2.7- 5.5	Accept.

Section No. 12 Revision No. 3 Date: 10/17/96 Page 11 of 18

Performance Evaluation Report USIPA Water Pollution Study WF035

rage: A Tate: 187886

Fartici	pant ID:	PAGCCO9		Type: OTHER	Request 1	na ((fice: U1
	Number	Reported Value	Value≄	Acceptance Limits	Warning Limits	Fraluation
I C B	'S TE CT	. IK PILLT	C03F2 \K1	100619		
		1016/1242		1 CC H H I		
,,, ,,,,	7. 012 C1	31.1	u2.3	6.00- 50.9	17.5- 67.7	topent.
01-PCH	IN CIL-		74.5			,
.or teb	02	14.9	12.7	3.17- 20.4	5.37- 10.2	iccept.
PES	TICICES	IN MICROGR	AMS/LITE	· F		
147-ALD			,			
	<b>C 1</b>	3.01	3.11	0.522- 5.23	1.12- 4.64	rccert.
		0.184	0.243	0.065-0.322	0.0977-0.789	Accept.
113-84C						,
	C1	0.38	4.51	2.67- 6.22	7.09- 5.76	Accept.
	(2	1.63	1.62	2.67- 6.22 0.058- 2.19	1.03- 2.02	Accept.
49-DDD		,			· · · -	
		6.10	5.67	1.14- 9.33	3.92- 8.55	Accent.
	(2	1.07	1.54	1.21- 2.64	1.79- 2.46	Accept.
50-DDE		100,		,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,		· · · · ·
,,,,,		1.60	1.76	2.14- 5.1	2.51- 0.73	*ccept.
	62	1.75	1.42	0.77- 1.85	0.063- 1.71	iccent.
51-DDT		2	••••			4.1,11
	(1	6.61	6.06	3.79- 9.28	11.40-0.59	*conpt.
	02	1.63	1.76	0.865- 2.33	1.05- 2.14	Accept.
152-88P	TACHLOR	200.,			20 20	
	C 1	2.41	2.85	C.694- 4.14	1.17- 2.71	Accept.
	02	0.232	0.278	C.694- 4.14 O.C899-0.374	C.126-C.138	Accept.
153-CBI	ORDANE	0				
,,,,		12.2	12 7	4.69- 17.2	6.27- 15.6	sacept.
	Č 4			0.695- 1.61	0.935- 1.67	*ccept.
78-HED	-	POTICE	1	01073 1101		acco per
,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,	01	1 80	2 20	1.17- 2.53	1 71- 2 75	Accent.
	(2	[.530	C - 280	0.153- 0.37	C.10-C.702	*ccent.
					et seture	£ 1, € 1, 3 €
	ATILE HAL		IA WICBO	GEAUSZLITER		
, J. 4	C1	7 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	56 7	37- 75.4	00.1- 70.6	tecert.
	02	13.4	12.2	P.5- 17.5	7.63- 16.4	Accept.
155_~#	OROFORM	3 .3 6 74	14.4	• •	2011 - 1501	acrapt.
ハンユイロF	C1	70.6	64.8	47.6- 83.2	12.1- 78.7	Accopt.
			14.2	11- 19.4	11.7- 17.4	*ccopt.
	02	13.8		) t= - t ( • 4	110 - 1100	icinpi.
7.1-0c		LOFOETHANE		41-06-7	H C C . OO 1	10000
	01	63.1	63.7	01- 05.7	46.6- 20.1	*ccept.
	C 2	17.1	16.2	10.4- 22.2	11.9- 20.7	Accept.
157-TRI	CHLORNET					
	CJ	7n.E	15.9	n : . 1 - 93	<1.5- N7.1.	topopt.
	( 2	16.3	16.1	10.7- 20.6	11.0- 19.4	≇ccert.
258-CAR	BONTETRA					
	C 1	79.7	39.C	16.5- 45	7C- 41.4	tecept.
	0.2	10.4	9.36	6.61- 13.7	6.67- 17.7	tecent.

Section No. 12 Revision No. 3 Date: 10/17/96 Page 12 of 18

# Performanco Evaluation Report USEPA Water Pollution Study WPC35

Page: 5 Date: 164ppg6

Participant ID:	PAOCCCO		Type: CIUEF	Pequest1	ng (ffice: UT
Ruster	Value	7cue Value⊅	Acceptance Iluits	]   m   t s	FTAINATICE
C59-TETRACII I ORO					
(1		73.6	46.5- 96.3	52.7- 99.1	Accept.
02	10.9	10.4	6.04- 14.5	7.1- 13.4	Accent.
C60-BROMOCICHLO					
	6C.2	55.6	37- 73.1	41.6- 6P.6	Accept.
	14.3	14.6			
061-DIB FORCCHLO	ROMETHANE				
01	51.0	48.5	33.1- 63.4	37- 59.6	locept.
(2	13.7	14.6	9.59- 19.9	10.7- 17.7	Accept.
162-BRONCFCFM					
		60.0	50.4- 95.3	56- PO.7	Accept.
	13.3	12.6	9.42- 17.3	9.53- 16.2	<pre></pre>
63-METHYLENE C	HLORIDE		30 B 44	3.11 6 6 6 6	<b>.</b>
C1			30.3- 64.1		
		10.3	6.63- 14.7	1.65- 11.1	Accept.
164-CHLCRCHENZE	NE TO A		n.r. c	E1 0 02 2	1,000
61	76.8	00.1	46.5- 89 11.7- 24.2	17 7 77 6	ACCEPT.
(2	17.8	1 / . /	31.7- 24.2	13.75 22.6	ACCRIT.
VCLATIIE AF 065-DENZENE					
(1	56.0	e ë " U	111.7- 69.5	44.4- 66.7	Accept.
0.2	9.25	9.30	6.56- 12.3	7.29- 11.6	iccort.
D66-ETHYLEFNZEN					
C 1	e7.4		36.7- 73.3	43- 69	Accert.
62	10.4	10.4	7.18- 13.5	7,99- 17,9	Accept.
167-1CLUENE					
61	44.6	44.7	30.9- 57.6	34.3- 54.2	
( 2	7.46	7.60	5.29- 9.97	5. No- 9.36	Accept.
)94-1,2-DICHLOR					
01	119.3		40.7- 66.5	44-63.2	
	11.0	11.7	7.02- 16.6	P. 92- 15.5	Accept.
)95-1,4-DICHLOR					
			33.3- 62.2		
		13.0	9.37- 17.6	10.4- 16.6	rccept.
96-1,3-DICHLOR			30 0 53 4	24 2 56 5	
01	37.1	42.7	30.4- 53.1	36.7- 50.7	tocept.
C 2	11.1	12.6	e.e- 16.6	9.70- 15.6	'ccept.
MISCELL/ NEO					
71-TCTAL CYANI	•	•			
01	.022	.0321	0.0139 - 0.046	0.0738-0.003	•
(2	•406	C.41C	0.297-0.522	0.352-0.443	Accept.
172-ACN-FILTERA		-			
91	64.6	86.0	(1.7- 28.1	66.4- 93.6	Ct. for Frr.
C 5	0(.0	56.0	17.7- 60.1	44.9- 57.9	Accept.
373-CIL AND GRE					
6.1	51.F	44.0	30 ° c = 20 ° T	72.6- 51	Charter for ter-
0.2	23.5	18.9	10- 23.2	13.4- 21.9	icanak.

Section No. 12 Revision No. 3 Date: 10/17/96 Page 13 of 18

Performance Evaluation Report USEPA Water Pollution Study WP035 Date: 16\*P\*06

rage: E

Farticipan	t ID:	PACCCCO		Type: OTHER	Request	ina Cffice: OT
S a N u	mple	Reported Value	True Value#	Acceptance Limits	Vaining Limits	Performance rvaluation
097-101AL	PHERCL	.ICS (IK MG	/f.)	,		
	C1			1.47- 3.96	1.78- 3.54	*crert.
	<b>C</b> 2	.011	1.19	0.519- 1.87	0.607- 1.7	Accept.
098-TCTAL	RESICU	IAL CFICFI	SECIN MO	5/1)		·
			•	7.50- 3.6	7.68- 7.46	Accept.
	( 2	.32C	C.41C	0.295-0.624	L.330-L.2E1	ct. for fer.
<b>,,,,,,,,</b> ,,,,,,,,,,,,,,,,,,,,,,,,,,,,	END C	OF DATA FO	R FACCC	9 ******		
NOTE: FOR	LINIT	S AND TRU	F VALUES	, ASSUNE THREE	SIGNIFICANT	DIGITS.
				0000 0000000000000000000000000000000000		

<sup>\*</sup> Based on gravimetric calculations, or a reference value when necessary.

Section No. 12 Revision No. 3 Date: 10/17/96 Page 14 of 18

#### ORGANIC PREAWARD EVALUATION SAMPLE INDIVIDUAL LABORATORY SUMMARY REPORT

LABORATORY: Lancaster Laboratories (PA) [LANCA2]

SCORE: 100.0

MATRIX: WATER DATE: 11/30/94

					DATE:	11/30/94	
		מספחוריים	N INTERVAL	c			
	GAW.	PREDICTION		JON	LABORATORY DATA		
COMPOUND	LOWER	UPPER	LOWER	UPPER	CONC	OUAL	
COMBOUND							
TCL VOLATILE							
VINYL CHLORIDE	43	92	35	120	74		
CHLOROFORM	59	77	56	79	72		
CARBON TETRACHLORIDE	76	100	72	110	100		
BROMODICHLOROMETHANE	63	79	61	81	80	2	
TRICHLOROETHENE	120	170	120	170	160		
DIBROMOCHLOROMETHANE	83	100	80	100	96		
BENZENE	77	<del>96</del>	75	98	84		
BROMOFORM	<b>87</b>	110	83	120	96		
TOLUENE	61	<b>7</b> 7	58	79	68		
XYLENES (TOTAL)	96	130	90	140	100		
TCL SEMIVOLATILE							
PHENOL	35	54	32	64	49		
4-METHYLPHENOL	24	40	21	49	28		
ISOPHORONE	24	35	23	36	30		
2.4-DIMETHYLPHENOL	12	37	0	50	. 9	S	
1,2,4-TRICHLOROBENZENE	25	39	23	47	37		
NAPHTHALENE	20	28	19	32	26		
HEXACHLOROCYCLOPENTADIENE	10	64	0	73	18		
2.4,6-TRICHLOROPHENOL	43	62	40	65	55		
4-BROMOPHENYL PHENYL ETHER	33	44	31	45	40		
HEXACHLOROBENZENE	27	36	26	37	33		
FLUORANTHENE	34	45	33	47	46	S	
PYRENE	40	59	37	62	44	•	
DI-N-OCTYL PHTHALATE	37	65	33	69	63		
BENZO(A)PYRENE	15	43	11	58	24		
DIBENZ(A,H)ANTHRACENE	34	58	30	62	44		
TCL PESTICIDES							
BETA-BHC	0.27	0.41	0.25	0.43	0.36		
HEPTACHLOR EPOXIDE	0.39	0,6	0,36	0.62	0.52		
DIELDRIN	0.42	0.65	0,38	0.68	0.54		
4.4'-DDE	0.36	0.61	0.32	0.64	0.51		
ENDRIN KETONE	0.6	1.1	0.53	1.1	0.87		
GAMMA-CHLORDANE	0.35	0.53	0.33	0.56	0.46		
		05	0.55	0,50	v. <del>4</del> 0		

Section No. 12 Revision No. 3 Date: 10/17/96 Page 15 of 18

#### ORGANIC PREAWARD EVALUATION SAMPLE INDIVIDUAL LABORATORY SUMMARY REPORT

LABORATORY: Lancaster Laboratories (PA) [LANCA2]

SCORE: 100.0

MATRIX: WATER

DATE: 11/30/94

#### PREDICTION INTERVALS

	WAR	NING	ACT	ION	LABORATO	RY DATA
COMPOUND	LOWER	UPPER	LOWER	UPPER	CONC.	OUAL
NON-TCL VOLATILE						
BENZENE.N-PROPYL- HEXANE PROFANE,1.2-DIBROMO-3-CHLORO-					63 62	NR
NON-TCL SEMIVOLATILE						
BENZYL ALCOHOL DIBENZOTHIOPHENE PARATHION					44 42 15	
TCL VOLATILE (Contaminants)						
METHYLENE CHLORIDE					2	
TCL SEMIVOLATILE (Contaminants)						
DIETHYLPHTHALATE					2	
TIC SEMIVOLATILE (Contaminants)						
UNKNOWN CHLORINATED COMPOUND					3	
# OF TCL COMPOUNDS NOT-IDENTIFIED: 0 # OF TCL COMPOUNDS MIS-QUANTIFIED: 0 # OF TCL CONTAMINANTS: 0						

<sup>#</sup> OF NON-TCL COMPOUNDS NOT-IDENTIFIED: 0

Section No. 12 Revision No. 3 Date: 10/17/96 Page 16 of 18

#### ORGANIC PREAWARD EVALUATION SAMPLE INDIVIDUAL LABORATORY SUMMARY REPORT

LABORATORY: Lancaster Laboratories (PA) [LANCA2]

**SCORE**: 93.3

MATRIX: SOIL DATE: 11/30/94

חשמ	TON	INTERVA	T	9

		<b>PREDICTIO</b>	N INTERVAL	S		
	WAR	NING	ACT	ION	LABORATO	RY DATA
COMPOUND	LOWER	UPPER	LOWER	UPPER	CONC	OUAL
TCL VOLATILE						
CHLOROETHANE	57	140	44	160	120	
1.2-DICHLOROETHANE	58	81	55	84	66	
1.2-DICHLOROPROPANE	56	80	53	84	72	
TRICHLOROETHENE	81	120	76	120	120	
DIBROMOCHLOROMETHANE	92	150	84	150	110	
BENZENE	74	110	69	110	94	
TETRACHLOROETHENE	76	100	72	110	88	
TOLUENE	53	<i>7</i> 3	50	76	55	
CHLOROBENZENE	100	140	<del>9</del> 7	140	130	
XYLENES (TOTAL)	94	130	88	140	90	S
TCL SEMIVOLATILE						
2-METHYLPHENOL	900	2000	740	2600	200	x
NITROBENZENE	970	2100	810	2600	1900	
2-METHYLNAPHTHALENE	1200	2300	1100	2900	1800	
2-CHLORONAPHTHALENE	1300	2500	1100	3200	2400	
ACENAPHTHYLENE	1600	2700	1400	3200	2300	
DIBENZOFURAN	1600	2800	1400	3500	2400	
FLUORENE	1500	2500	1400	3000	1800	
PENTACHLOROPHENOL	3000	7600	2300	10000	5900	
PHENANTHRENE	1.700	3100	1500	3800	2500	
ANTHRACENE	1700	2800	1500	3400	2100	
CARBAZOLE	1400	2800	1200	3000	2200	
CHRYSENE	1600	2700	1400	2800	2000	
BENZO(B)FLUORANTHENE	1300	2300	1200	2400	2000	
INDENO(1,2,3-CD)PYRENE	1200	2200	1100	2700	1600	
BENZO(G.H,I)PERYLENE	1200	2900	920	3800	1900	
TCL PESTICIDES						
GAMMA-BHC (LINDANE)	8.7	23	6.6	25	18	
HEPTACHLOR	12	26	9.6	28	21	
ALDRIN	6.7	17	5.3	18	12	
ENDRIN	27	59	22	63	46	
ENDOSULFAN II	19	49	15	53	33	
4,4'-DDT	20	51	15	56	38	
METHOXYCHLOR	120	280	96	300	250	

Section No. 12 Revision No. 3 Date: 10/17/96 Page 17 of 18

#### ORGANIC PREAWARD EVALUATION SAMPLE INDIVIDUAL LABORATORY SUMMARY REPORT

LABORATORY: Lancaster Laboratories (PA) [LANCA2]

SCORE: 93.3 MATRIX: SOIL DATE: 11/30/94

PREDI	CTI	ON:	ME	RVALS
-------	-----	-----	----	-------

COMPOUND	WARNING		ACTION		LABORATORY DATA	
	LOWER	UPPER	LOWER	UPPER		QUAL.
CON EQ.						
NON-TCL VOLATILE						
BENZENE,N-PROPYL-					50	
HEXANE					72	
PROPANE, 1.2-DIBROMO-3-CHILORO-						NR
NON-TCL SEMIVOLATILE						
DELITE DA					5300	
BIPHENYL DIPHENYL HYDRAZINE					2000	
PARATHION					1600	
FARATHION					1000	
TCI. VOLATILE (Contaminants)						
METHYLENE CHLORIDE					3	
ACETONE					52	
2-BUTANONE					10	
TCI. SEMIVOLATILE (Contaminants)						
ACENAPHTHENE					100	
BIS(2-ETHYLHEXYL)PHTHALATE					310	
TCL PESTICIDES (Contaminants)						
DEL.TA-BHC					0.33	
ENDOSULFAN I					0.26	
4.4'-DDE			-		0.33	
ENDRIN KETONE					2.3	
ENDRIN ALDEHYDE					1.4	
TIC VOLATILE (Contaminants)						
ETHANE,1,1,2-TRICHLORO-1,2,					11	
TIC SEMIVOLATILE (Contaminants)						
9,10-ANTHRACENEDIONE					=0	
UNKNOWN					78 <b>8</b> 2	
UNENOWN					82 75	
wertenser with the					13	

Section No. 12 Revision No. 3 Date: 10/17/96 Page 18 of 18

### ORGANIC PREAWARD EVALUATION SAMPLE INDIVIDUAL LABORATORY SUMMARY REPORT

LABORATORY: Lancaster Laboratories (PA) [LANCA2]

SCORE: 93.3

MATRIX: SOIL

DATE: 11/30/94

PREDICTION INTERVALS

WARNING ACTION LABORATORY DATA
COMPOUND LOWER UPPER LOWER UPPER CONC. QUAL

# OF NON-TCL COMPOUNDS NOT-IDENTIFIED: 0

<sup>#</sup> OF TCL COMPOUNDS NOT-IDENTIFIED: 0

<sup>#</sup> OF TCL COMPOUNDS MIS-QUANTIFIED: 1

<sup>#</sup> OF TCL CONTAMINANTS: 0

Section No. 13 Revision No. 2 Date: 07/25/95 Page 1 of 2

### 13. Preventive Maintenance

In order to ensure timely production of data, Lancaster Laboratories schedules routine preventive maintenance of instruments based on manufacturer's recommendations. Maintenance of the laboratory instruments is the responsibility of the technical group using the equipment in conjunction with our in-house Equipment Maintenance Group. A schedule of routinely performed instrument rnaintenance tasks is attached as Table 13-1. All preventive maintenance, as well as maintenance performed as corrective action, is recorded in instrument logs.

Critical spare parts are kept in supply at the laboratory by the Equipment Maintenance Group. Most items not kept in stock at the laboratory are available through overnight delivery from the manufacturer. In addition, Lancaster Labs maintains multiple numbers of most of the critical instruments used in our laboratory operations. A recent equipment inventory may be found in the *Qualification Manual*. Because we are a large laboratory with redundant capacity, the problems of instrument downtime are minimized.

Section No. 13 Revision No. 2 Date: 07/25/95 Page 2 of 2

	Table 13-1	
	Preventive Maintenance Schedule	
Instrument	Preventive Maintenance Schedule Preventive Maintenance	Frequency
GC/MS	Column maintenance	AN*
40/11/0	Change septum	Weekly or AN*
	Check fans	Monthly
	Check cool flow	Monthly
	Clean source	Bimonthly or AN
	Change oil in vacuum pump	Semiannually
	Change oil in turbo pump	Semiannually
GC	Septum change	Each run
40	Column maintenance	AN
	Clean detector	AN
	Vacuum filters	Semiannually
	Leak check ECDs	Semiannually
Flame AA	Rinse burner head, chamber and trap	AN: Min. Weekly
	Clean nebulizer	Weekly
	Inspect tubing and O-rings	Monthly
	Replace lamp	AN
GFAA	Rinse workhead assembly	Weekly
	Clean windows	Weekly
	Replace probe tubing	AN
	Check rinse bottle & drain	Daily
ICP	Clean torch	AN
	Clean nebulizer & spray chamber	AN
	Replace pump winding	Check Daily
	Lubricate autosampler	Check Daily
	Check mirror	Check Daily
	Check tubing to torch	Daily
	Check fan filters, clean if needed	Weekty
	Check cool flow, clean if needed	Weekly
	Check water filter, replace if needed	Quarterly
Cold Vapor AA	Change drying tube	Daily
•	Replace pump tubing	AN: Min. weekly
	Lubricate pump head	Weekly
	Lubricate autosampler	Weekly
	Inspect optical cell and windows	Monthly
	Clean	AN
Autoanalyzer	Clean sample probe	AN
•	Clean proportioning pump	Weekly
	Inspect pump tubing, replace if worn	AN
	Clean wash receptacles	Monthly
<u> </u>	Inspect condition of distillation head	Monthly

<sup>\*</sup>AN means as needed. Any of these items may be performed more frequently if response during operation indicates this is necessary.

Section No. 14 Revision No. 2 Date: 07/25/95 Page 1 of 3

# 14. Specific Routine Procedures Used to Assess Data Precision, Accuracy, and Completeness

<u>Precision</u> - Precision refers to the reproducibility of a method when it is repeated on a second aliquot of the same sample. The degree of agreement is expressed as the relative percent difference (RPD). The RPD will be calculated according to the following equation:

$$RPD = \frac{D_2 - D_1}{(D_1 + D_2)/2} \times 100$$

Where:

 $D_1$  = First sample value

 $D_2$  = Second sample value (Duplicate)

Duplicates will be run on at least 5% of the samples. Acceptance criteria shall be within the value range specified by EPA in the CLP SOW. (See Section No. 11.) All quality control sample results are entered into the computer and compared with acceptance limits. In addition, there is a monthly review of values on the computer QC system. Data obtained from quality control samples is entered onto our computer system which charts the data and calculates a mean and standard deviation on a monthly basis. The Quality Assurance Department then reviews this data for trends which may indicate analytical problems. The control charts are graphical methods for monitoring precision and bias over time.

<u>Accuracy</u> - Accuracy refers to the agreement between the amount of a compound measured by the test method and the amount actually present. Accuracy is usually expressed as a percent recovery (R). Recoveries will be calculated according to the following equations:

Section No. 14 Revision No. 2 Date: 07/25/95 Page 2 of 3

Surrogate Re cov ery = 
$$\frac{Qd}{Qa} \times 100$$

Where:

Qd = Quantity determined by analysis

Qa = Quantity added to sample

Matrix Spike Re cov ery = 
$$\frac{SSR - SR}{SA} \times 100$$

Where:

SSR = Spiked sample results

SR = Sample results

SA = Spike added

Laboratory Control Sample Re cov ery = 
$$\frac{LCS Found}{LCS True} \times 100$$

Surrogate standards are added to each sample analyzed for organics. Spikes and laboratory control samples will be run on at least 5% of the samples (each batch or SDG, ≤20 samples). Acceptance criteria for the accuracy recoveries shall be within the range specified by EPA in the CLP SOW. (See Section No. 11.) The Laboratory computer is programmed to compare the individual values with the acceptance limits and inform the analyst if the results meet specification. If the results are not within the acceptance criteria, corrective action suitable to the situation will be taken. This may include, but is not limited to, checking calculations and instrument performance, reanalysis of the associated samples, examining other QC analyzed with the same batch of samples, and qualifying results with documentation of any QC problems in the case narrative.

Commercial quality control materials are run at least quarterly to ensure accuracy of the analytical procedure. Repetitive analysis of a reference material will also yield precision data. Accuracy information determined from reference materials is valuable because variables specific to sample matrix are eliminated.

Section No. 14 Revision No. 2 Date: 07/25/95 Page 3 of 3

The QC program is capable of charting data for surrogates, spikes, control materials, and reference materials. The Quality Assurance Department reviews these charts for any indication of possible problems (i.e., shift in the mean and standard deviation).

Completeness - Completeness is the percentage of valid data acquired from a measurement system compared to the amount of valid measurements that were planned to be collected. The objective is analysis of all samples submitted intact, and to ensure that sufficient sample weight/volume is available should the initial analysis not meet acceptance criteria. The laboratory's sample management system will assign a unique identification number to the sample which tracks and controls movement of samples from the time of receipt until disposal. All data generated will be recorded referencing the corresponding sample identification number. The completeness of an analysis can be documented by including in the data deliverables sufficient information to allow the data user to assess the quality of the results. This information will include, but is not limited to, summaries of QC data and sample results, chromatograms, spectra, and instrument tune and calibration data. Additional information will be stored in the laboratory's archives, both hard copy and magnetic tape.

 $Completeness = \frac{Number of valid measurements}{Total measurements needed} \times 100$ 

Section No. 15 Revision No. 2 Date: 07/25/95 Page 1 of 3

### 15. Corrective Action

Whenever any of the data generated falls outside of the established acceptance criteria outlined for instrument tune and calibration (Section 8) and internal QC (Section 11), the cause of this irregularity must be investigated, corrected, and documented. The documentation will be used to prevent a recurrence of the problem and to inform management of the situation.

If the results are not within acceptance criteria, the appropriate corrective action will be initiated. This may include, but is not limited to, checking calculation and instrument performance, reanalysis of the associated samples, examining other QC analyzed with the same batch of samples, and qualifying results with a comment stating the observed deviation.

A standard operating procedure is in place which outlines the procedures to be followed when quality control data for an analysis falls outside of previously established acceptance limits. All QC data must be entered onto the computerized QC system promptly after its generation and daily "out-of-spec" data is reported via this system. Any data outside the acceptance criteria will be reviewed by the Quality Assurance Department. Where appropriate, the Quality Assurance Department will place outliers in one of three categories:

- A. <u>Marginal Outlier</u> Data that are outside the 95% confidence interval but within the 99% confidence interval. This category may also be used for QC samples subject to matrix interferences or sample inhomogeneity.
- B. Outlier Data outside the 99% confidence interval and/or observable trends such as a shift in mean and standard deviation.
- C. <u>Extreme Outlier</u> Such data would indicate the system is out of control and no results should be reported to clients; an example would be more than one reference or control falling outside the 99% confidence interval.

The daily out-of-spec reports are then distributed to group leaders or their QC coordinator who will check all supporting data and document their findings and any corrective action taken. Documentation of QC data will be filed in the departmental QC notebook. In the case of outliers or extreme outliers, the Quality

Section No. 15 Revision No. 2 Date: 07/25/95 Page 2 of 3

Assurance Department may issue a formal request for investigation and corrective action (see sample form that follows). The Quality Assurance Department is responsible for initiating the corrective actions, insuring that the actions are taken in a timely manner, and that the desired results are produced. The QA Department will circulate all completed Investigation & Corrective Action forms to the appropriate manager.

The Quality Assurance Department is also responsible for conducting periodic audits which ensure compliance with laboratory SOPs and assist in identifying and correcting any deficiencies. These audits may entail observation as procedures are carried out or a review of records to demonstrate traceability and compliance with all documented record keeping procedures. The QA Department will issue a audit written report which summarizes the audit findings and the technical centers are then requested to respond in writing within 30 days of report receipt. The response will address the corrective action that needs to be taken along with an expected completion date. Audit results and the corresponding response are communicated to laboratory personnel and management. Follow-up audits verify that proper corrective action has been taken for the identified discrepancy.

Section No. 15 Revision No. 2 Date: 07/25/95 Page 3 of 3



No.		

Investigation and Corrective Action Report
Part I Description of problem
<ol> <li>Date</li> <li>LLI sample number(s) involved</li> <li>Nature of problem (e.g., QA outlier, procedural deviation, client complaint, etc.)</li> </ol>
4. Check if investigation must be complete before reporting further data to clients
Initiated by:
Part II (Attach separate sheet if needed)
1. Steps taken to investigate problem.
2. Explanation of probable cause of problem.
3. Steps taken to prevent future occurrence.
4. Besides the sample(s) listed above, would data sent to any clients be affected by this problem? If yes, explain.
5. Signed: Date:
Return by:

2064 Rev. 09/01/94

Section No. 16 Revision No. 2 Date: 07/25/95 Page 1 of 1

### 16. Quality Assurance Reports to Management

Reports of quality status from the Quality Assurance Department to management are made frequently and in various forms. All results from internal or external performance evaluation samples are circulated to management. A report of each audit performed is prepared and copied to management. Monthly summaries of data obtained from analysis of quality control check samples are generated via the computerized sample management system. These summaries include mean and standard deviation to aid in assessment of data accuracy and precision. Forms summarizing problems which require investigation and corrective action are completed by group leaders and circulated to management. Through these channels, laboratory management is kept apprised of QA/QC activities.

Any problems or unusual observations that occur during the analysis of samples for a specific project will be listed on the laboratory report and/or in the case narrative delivered with the data package. The items often discussed in this manner include samples with surrogate recovery outside of the acceptance criteria and samples with matrix problems requiring dilution and causing increased detection limits. Where applicable, any corrective action attempted or performed to address the problem will also be presented.

The laboratory will contact the client for direction regarding major problems such as samples listed on the chain of custody but missing from the shipping container, samples which arrive broken or are accidentally broken in the laboratory, and samples with severe matrix problems. The client will be contacted if it is necessary to change any item in the original project plan.

Appendix A

CLP Forms
Inorganics and Organics

### VOLATILE ORGANICS ANALYSIS DATA SHEET

Lab Name: Co	ontract:	
Lab Code: Case No.: S	AS No.:SDG No.:	
Matrix: (soil/water)	Lab Sample ID:	
Sample wt/vol:(g/mL)	Lab File ID:	_
Level: (low/med)	Date Received:	
% Moisture: not dec	Date Analyzed:	
GC Column: ID:(mm)	Dilution Factor:	
Soil Extract Volume:(uL)	Soil Aliquot Volume:(uL	( د
CAS NO COMPOUND	CONCENTRATION UNITS: (ug/L or ug/Kg)Q	
74-87-3	de  ne ne ne ne (total)  ne thane ride hane ane propene hane thane thane	
100-41-4Ethylbenzene 100-42-5Styrene 1330-20-7Xylene (total)		
	1	

### SEMIVOLATILE ORGANICS ANALYSIS DATA SHEET

		•	. 1	
.ab Name:		Contract:		
ab Code:	Case No.:	SAS No.:	SDG No.	. :
Matrix: (soil/wate	r)	Lab	Sample ID: _	
Sample wt/vol:	(g/mL)	Lab	File ID:	*_**_*********************************
evel: (low/med)		Date	Received:	
Moisture:	decanted: (Y/N)_	Date	Extracted:_	<del></del>
oncentrated Extra	ct Volume:(	uL) Date	Analyzed:	
njection Volume:	(uL)	Dilu	tion Factor:	
PC Cleanup: (Y/	N) pH:	•		
CAS NO.	СОМРОИИД		CION UNITS:	Q
108-95-2	Phenol			
	bis(2-Chloroet	hvl)ether	<b>-</b>	
95-57-8	2-Chlorophenol			<u> </u>
541-73-1	1,3-Dichlorobe	nzene	-	<del>-</del>
106-46-7	1,4-Dichlorobe			<del>-</del>
95-50-1	1,2-Dichlorobe	nzene		<del></del>
1 95-48-7	2-Methvlnhenol		1	-
108-60-1	2,2'-oxybis(1-	Chloropropane	5	— <u>  —                                  </u>
106-44-5	4-Methylphenol		<b>'</b>	<u> </u>
621-64-7	N-Nitroso-di-n	-propylamine	-	<del></del>
67-72-1	Hexachloroetha	ne		_
98-95-3	Nitrobenzene		1	1 1
78-59-1	Isophorone	<del></del>	_	-
88-75-5	2-Nitrophenol		4	
105-67-9	2,4-Dimethylph	enol	-	<del></del>
111-91-1	bis(2-Chloroet	hoxy)methane	-	
120-83-2	2,4-Dichloroph	enol	<del>-</del>	_
120-82-1	1,2,4-Trichlor	obenzene	_	
	Naphthalene			_
106-47-8	4-Chloroanilin			
	Hexachlorobuta			
59-50-7	4-Chloro-3-met	hylphenol		
91-57-6	2-Methylnaphth	alene		
77-47-4	Hexachlorocycle	opentadiene		
88-06-2	2,4,6-Trichlor	ophenol		
95-95-4	2,4,5-Trichlore	ophenol		
91-58-7	2-Chloronaphth	alene		_
	2-Nitroaniline		-	
	Dimethylphthal		*	_
	Acenaphthylene			_
	2,6-Dinitrotol		_	_
	3-Nitroaniline		-	— <u> </u> ——-
	Acenaphthene		~	_
			-	— ——

#### 1C SEMIVOLATILE ORGANICS ANALYSIS DATA SHEET

Lab Name:		Contract:		
Lab Code:	Case No.:	SAS No.:	SDG No	
Matrix: (soil/water)	)	Lab	Sample ID: _	
Sample wt/vol:	(g/mL)	Lab	File ID: _	<del></del>
Level: (low/med)		Date	Received: _	<del></del>
% Moisture:	decanted: (Y/N	) Date	Extracted:_	
Concentrated Extract	: Volume:	_(uL) Date	Analyzed: _	
Injection Volume:	(uL)	Dilu	tion Factor:	
GPC Cleanup: (Y/N)	pH:	<del>-</del>	•	•
CAS NO.	СОМРОИЙД	CONCENTRAT	ION UNITS: g/Kg)	Q
100-02-7 132-64-9 121-14-2 84-66-2 7005-72-3 86-73-7 100-01-6 534-52-1 86-30-6 101-55-3 18-74-1 87-86-5 85-01-8 120-12-7 86-74-8 206-44-0 129-00-0 85-68-7 91-94-1 56-55-3 218-01-9 117-84-0 205-99-2 207-08-9 50-32-8 193-39-5	Butylbenzylph 3,3'-Dichlord Benzo(a)anthr	cluene late yl-phenylether ne 2-methylphenol nenylamine (1) l-phenylether nzene nenol chalate benzidine racene exyl)phthalate chalate ranthene ne ned)pyrene nethracene		

(1) - Cannot be separated from Diphenylamine

. 1D

EPA SAMPLE NO.

### PESTICIDE ORGANICS ANALYSIS DATA SHEET

		<b>-</b>		i v
Lab Name:		Contract:		
Lab Code:	Case No.:	SAS No.:	SDG No.:	<del></del>
Matrix: (soil/water)		Lab S	ample ID:	
Sample wt/vol:	(g/mL)	Lab F	ile ID:	
Moisture:	decanted: (Y/N)_	Date	Received:	
Extraction: (SepF/C	ont/Sonc)	_ Date	Extracted:	·
Concentrated Extract	Volume:(	uL) Date	Analyzed:	<del></del>
Injection Volume:	(uL)	Dilut	ion Factor: _	<del></del>
GPC Cleanup: (Y/N)	pH:	Sulfu	r Cleanup: (Y,	/N)
CAS NO.	COMPOUND	CONCENTRATI		Q
319-85-7 319-86-8 58-89-9 76-44-8 309-00-2 1024-57-3 959-98-8 72-55-9 72-20-8 33213-65-9 72-54-8 1031-07-8 50-29-3 72-43-5 53494-70-5 7421-93-4 5103-71-9 5103-74-2 8001-35-2 12674-11-2 11104-28-2 11104-28-2 11097-69-1	Endosulfan II4,4'-DDDEndosulfan sul4,4'-DDTMethoxychlorEndrin ketoneEndrin aldehyddalpha-Chlordand	dane) xide fatee		

# VOLATILE ORGANICS ANALYSIS DATA SHEET TENTATIVELY IDENTIFIED COMPOUNDS

Lab Name:	Contract:
Lab Code: Case No.:	SAS No.:SDG No.:
Matrix: (soil/water)	Lab Sample ID:
Sample wt/vol:(g/mL)	Lab File ID:
Level: (low/med)	Date Received:
% Moisture: not dec	Date Analyzed:
GC Column: ID: (mm)	Dilution Factor:
Soil Extract Volume:(uL)	Soil Aliquot Volume:(uL)
Number TICs found:	CONCENTRATION UNITS:

CAS NUMBER	COMPOUND NAME	RT	EST. CONC.	
1				
2				
3				
4				
5				
6	<del></del>			
7		.	·	
8				_
9				_
.0		-		_
· + •  _		.		
.2.		·		
. * •		.		
. 3 • [_		-		
.6				
	<del></del>	:		
.o		·		
		·		
0	<del></del>	·		—
1				<del></del>
2		·		
3	<del></del>			
4				
5				
6				
7				
8				
9				

1F

EPA SAMPLE NO.

# SEMIVOLATILE ORGANICS ANALYSIS DATA SHEET TENTATIVELY IDENTIFIED COMPOUNDS

Lab Name:	Contract:
Lab Code: Case No.:	SAS No.: SDG No.:
Matrix: (soil/water)	Lab Sample ID:
Sample wt/vol:(g/mL)	Lab File ID:
Level: (low/med)	Date Received:
% Moisture: decanted: (Y/N)	Date Extracted:
Concentrated Extract Volume:(	L) Date Analyzed:
Injection Volume:(uL)	Dilution Factor:
GPC Cleanup: (Y/N) pH:	
Number TICs found:	CONCENTRATION UNITS: (ug/L or ug/Kg)

1			EST. CONC.	Q
1.			=======================================	====
		1		
2				
3.		-		
3	· · · · · · · · · · · · · · · · · · ·		<del></del>	
4		·		
5				<u> </u>
6				
/·				
8				
9				
	<del></del>			
0				
• 📤 •	<del> </del>		i <del></del> i	
.4 •	······································	.		
		1	l	
4				
.5				
	<del></del>			
6	<del></del>			
. / •				
			<u> </u>	
• •				
0.				
1				
1	······································	·		
2				
J.				
4.		l		
6				
7				
· · ·				
9				

2A WATER VOLATILE SYSTEM MONITORING COMPOUND RECOVERY

Lab	Name:		Contract:	<del></del>
Lab	Code:	Case No.:	SAS No.:	SDG No.:

1	EPA	SMC1	SMC2	SMC3	OTHER	TOT
ł	SAMPLE NO.	(TOL)#	(BFB)#	(DCE)#	Jornan .	OUT
ŀ	STATE NO.	(102)		(505)#		===
01		-				
02	<del></del>	-				<b>—</b>
03		-				_
04	<del></del>	-	<del></del>		<del></del>	
05		-				
06		-				
07		-			ļ	
08		-				<b> </b> —
09	<del></del>	-			<b></b>	—
	<del></del> -	-		<del></del>		
10		_				<b>-</b>
		-				
12		-				
13		-				
14		-				_
15	·	-1				—
16		-				
17		-				
18	<del></del>	-				
19	<del></del>					
20	<del></del>				<del></del>	
21	<del></del>	_				
22		-		;		
23	· · · · · · · · · · · · · · · · · · ·	_				
24		-			·	
25	<del></del>	_				
26		_				
27						
28						
29		_				
30		_				

SMC1 (TOL) = Toluene-d8 (88-110) SMC2 (BFB) = Bromofluorobenzene (86-115) SMC3 (DCE) = 1,2-Dichloroethane-d4 (76-114)

# Column to be used to flag recovery values

\* Values outside of contract required QC limits

page \_\_ of \_\_

# 2B SOIL VOLATILE SYSTEM MONITORING COMPOUND RECOVERY

ab Name:			Contract	t:	<del></del>	
b Code:	_ Case No.:	<del></del>	SAS No.	• :	SD	G No.:
vel:(low/med)	·					
	EPA	SMC1	SMC2	SMC3	OTHER	TOT
i	SAMPLE NO.	(TOL)#	(BFB)#	(DCE)#	****	OUT
01						
02						
03						
04	<del></del>					
. OS . O6		<b></b>		<b></b>	ļ	
07	<del></del>					
08				<del></del>		
09						
10						
11						
12   13					<u> </u>	
14						
15						
16						
17						
18.						
19 20						
21						
22						
23						
24						
25	<del></del>				<del></del>	
26 27						
28						
29						
30						
ate to at	11 (MOT) = M-1	مور مسمدا	•		C LIMI	
	:1 (TOL) = Tol :2 (BFB) = Bro				(84-138) (59-113)	
	(BCE) = 1,2				70-121	

page \_ of \_

\* Values outside of contract required QC limits

N	ame:			C	ontract	:				
C	od <b>e:</b>	_ Case	No.: _	<del></del>	SAS No.	:	SDG	No.: _	<del></del>	
, •		r=		1 <del></del>	1 -4	\			T	7
	EPA SAMPLE NO.			' S3		S5 (2FP)#			S8 (DCB)#	
	*******	****	****	======	****	=====	=====	****	=====	
١.	· · · · · · · · · · · · · · · · · · ·				<u> </u>		<u> </u>	<u> </u>		
2].							<b> </b>			
3						<u> </u>				
5   ·	<del></del>			<del></del>						,
6 .	···		<del></del>	<del></del>					<del></del>	•
7	·									,
8				·						
9	<del></del>		<u> </u>							
9 .										,
1. 2								ļ		,
3										
4		<del></del>								
5   _										
6										,
7										
3 -										
3 -										
-										
2										
3										
1 _			<u> </u>							
5   -										
5 _										i
7   - 3   -										
<b>-</b>										
,   -										
							LIMITS			
		(NBZ)					5-114)			
		(FBP)	= Z-Fiu				3-116) 3-141)			

S4 (PHL) = Phenol-d5 (10-110)S5 (2FP) = 2-Fluorophenol (21-110)S6 (TBP) = 2,4,6-Tribromophenol (10-123) (33-110) (advisory) S7 (2CP) = 2-Chlorophenol-d4 S8 (DCB) = 1,2-Dichlorobenzene-d4 (16-110) (advisory) # Column to be used to flag recovery values

- \* Values outside of contract required QC limits
- D Surrogate diluted out

page \_ of \_

#### 2D SOIL SEMIVOLATILE SURROGATE RECOVERY

					·			
de:	_ Case	No.: _	***************************************	SAS NO.		SDG	мо.: _	
(low/med)								
·								
EPA	S1	S2	<b>S</b> 3	<b>S4</b>	S5	56	57	s
					(2FP)#		1	l .
######################################	22222	32222	#### <b>#</b>	*====	322222	*****		######
			<u> </u>					ļ
	<del></del>							
							] <del></del>	
							<del></del>	
	<u>-</u> _							
<del></del>								<b> </b>
					<del></del>			
			·					· ·
								ļ
								- <del></del>
<del></del>								
		<del></del>						
·						LIMITS		-

S1 (NBZ) = Nitrobenzene-d5 (23-120)

S2 (FBP) = 2-Fluorobiphenyl (30-115)

S3 (TPH) = Terphenyl-d14 (18-137)

S4 (PHL) = Phenol-d5 (24-113)

S5 (2FP) = 2-Fluorophenol (25-121)

S6 (TBP) = 2,4,6-Tribromophenol (19-122)

S7 (2CP) = 2-Chlorophenol-d4 (20-130) (advisory)

S8 (DCB) = 1,2-Dichlorobenzene-d4 (20-130) (advisory)

# Column to be used to flag recovery values

- \* Values outside of contract required QC limits
- D Surrogate diluted out

page \_ of \_

# 2E WATER PESTICIDE SURROGATE RECOVERY

ıb Code:		Case No.	· :	SAS	No.: _		SDG No	• :
: Column/	(1):	ID:	•	(mm) G	C Column	n(2):		rn:
CO Lumin,			·	(11411)	J 0014			
ŀ	EPA	TCX 1	TCX 2	DCB 1	DCB 2	OTHER	OTHER	TOT
	SAMPLE NO.	REC #	REC #	REC #	REC #	(1)	(2)	OUT
1	*******	*****	22222	****		22222		===
01								
02				ļ			<b> </b>	
03 04	<del></del>				<b></b>	<del></del>	<del></del>	
05						<del></del>		<del> </del>
06								
07						<del></del>		
08								
09								
10		•						
11								
12	<del></del>							
13								
. 14	<del></del>		·					
15								
16 17			<del></del>				<del></del>	
18								
19	<del></del>							<del></del>
20								
21	<del></del>				<del></del>			<del></del>
22								
23								
24								
25								
26								
27								
28								
29 30								

QC LIMITS

TCX = Tetrachloro-m-xylene (30-150) DCB = Decachlorobiphenyl (30-150)

# Column to be used to flag recovery values

\* Values outside of QC limits

D Surrogate diluted out

page \_\_ of \_\_

## SOIL PESTICIDE SURROGATE RECOVERY

Lab Name:		·		Cont	ract:				
Lab Code:		Case No	·:	SAS	No.: _	<del>,</del>	SDG No	.:	
GC Column	(1):	ID	·	(mm) G	C Columi	n(2): _	· · · · · · · · · · · · · · · · · · ·	ID:	(mm)
	EPA						OTHER	TOT	
	SAMPLE NO.	REC #	%REC #	%REC #	%REC #	(1)	(2)	OUT	•
	*******	32222	*****		250222	*****	****	===	
01					]				
02									
03									
04									
05									
06							<del></del>		
07									
08									
09									
10									
11									•
12									
13									
14									
15									
. 16							·		
17									
18									
19			·						
20									
21								<del></del>	
22								<del></del>	
23									
24									
4								<b></b> [	
25									
26									
27									
28									
29									
. 30									

### QC LIMITS

TCX = Tetrachloro-m-xylene (30-150)DCB = Decachlorobiphenyl (30-150)

# Column to be used to flag recovery values
\* Values outside of QC limits

D Surrogate diluted out

page or	page	of	
---------	------	----	--

#### 3A

### WATER VOLATILE MATRIX SPIKE/MATRIX SPIKE DUPLICATE RECOVERY

	CDIVE	CAVATA	- ve		
	SPIKE ADDED	SAMPLE CONCENTRATION	MS CONCÉNTERTOR	MS &	QC.
COMPOUND	(ug/L)	(ug/L)	(ug/L)		REC.
		( ug/ b/	( ug/ b)	RECF	AEC.
1,1-Dichloroethene					61-14
Trichloroethene				·	71-12
Benzene				·   <del></del>	76-12
Toluene					76-12
Chlorobenzene					75-13
COMPOUND  1,1-Dichloroethene	SPIKE ADDED (ug/L)	MSD CONCENTRATION (ug/L)	MSD REC # RPD #	RPD	IMITS   REC.
1,1-Dichloroethene Trichloroethene	ADDED (ug/L)	CONCENTRATION (ug/L)	% % REC ≠ RPD ≠	RPD	REC. ===== 61-14 71-12
1,1-Dichloroethene	ADDED (ug/L)	CONCENTRATION (ug/L)	% % REC ≠ RPD ≠	RPD 14 14 11	REC. ===== 61-14 71-12 76-12
	ADDED (ug/L)	CONCENTRATION (ug/L)	% % REC ≠ RPD ≠	RPD ====== 14 14	REC. ===== 61-14 71-12

#### 3*B*

### SOIL VOLATILE MATRIX SPIKE/MATRIX SPIKE DUPLICATE RECOVERY

.1		Contract: .				
ab Code: Case 1	No.:	_ SAS No.:	S	DG No.	:	_
atrix Spike - EPA Sample	No.:	Lev	el:(low/	med) _	<del></del>	
				_	·	
<del></del>	SPIKE	SAMPLE	MS		MS	QC.
:	ADDED	CONCENTRATION	CONCENT	RATION	*	LIMIT
COMPOUND	(ug/Kg)	(ug/Kg)	(ug/	Kg)	REC #	REC.
**********		*********	*****	******		1
1,1-Dichloroethene						59-17
Trichloroethene						62-13
Benzene		ļ				66-14
Toluene						59-13
Chlorobenzene						60-13
<u>.</u>						
	SPIKE	MSD	MSD			
	SPIKE ADDED	MSD CONCENTRATION	1 1	8	QC L	IMITS
COMPOUND			1 1	% RPD #	_	
*************	ADDED (ug/Kg)	CONCENTRATION	% 'REC ≠	RPD #	RPD	REC.
1,1-Dichloroethene	ADDED (ug/Kg)	CONCENTRATION (ug/Kg)	% 'REC ≠	RPD #	RPD ====== 22	REC.
1,1-Dichloroethene	ADDED (ug/Kg)	CONCENTRATION (ug/Kg)	% 'REC ≠	RPD #	RPD ====== 22 24	REC. 59-17 62-13
1,1-Dichloroethene Trichloroethene Benzene	ADDED (ug/Kg)	CONCENTRATION (ug/Kg)	% 'REC ≠	RPD #	RPD 22 24 21	REC. 59-17 62-13 66-14
1,1-Dichloroethene Trichloroethene Benzene Toluene	ADDED (ug/Kg)	CONCENTRATION (ug/Kg)	% 'REC ≠	RPD #	RPD 22 24 21 21	REC. 59-17 62-13 66-14 59-13
1,1-Dichloroethene Trichloroethene Benzene	ADDED (ug/Kg)	CONCENTRATION (ug/Kg)	% 'REC ≠	RPD #	RPD 22 24 21	REC. 59-17 62-13 66-14 59-13
1,1-Dichloroethene Trichloroethene Benzene Toluene	ADDED (ug/Kg)	CONCENTRATION (ug/Kg)	% 'REC ≠	RPD #	RPD 22 24 21 21	REC. 59-17 62-13 66-14 59-13
1,1-Dichloroethene Trichloroethene Benzene Toluene	ADDED (ug/Kg)	CONCENTRATION (ug/Kg)	% 'REC ≠	RPD #	RPD 22 24 21 21	REC. 59-17 62-13 66-14 59-13
1,1-Dichloroethene Trichloroethene Benzene Toluene Chlorobenzene	ADDED (ug/Kg)	CONCENTRATION (ug/Kg)	'REC #	RPD #	RPD 22 24 21 21 21	REC. 59-17 62-13 66-14 59-13
1,1-Dichloroethene Trichloroethene Benzene Toluene Chlorobenzene	ADDED (ug/Kg)	CONCENTRATION (ug/Kg)	'REC #	RPD #	RPD 22 24 21 21 21	REC. 59-17 62-13 66-14 59-13
1,1-Dichloroethene Trichloroethene Benzene Toluene Chlorobenzene  Column to be used to fla	ADDED (ug/Kg)	CONCENTRATION (ug/Kg)	'REC #	RPD #	RPD 22 24 21 21 21	REC. 59-17 62-13 66-14 59-13
1,1-Dichloroethene Trichloroethene Benzene Toluene	ADDED (ug/Kg)  servers ag recovery	CONCENTRATION (ug/Kg)  ===================================	REC #	RPD #	RPD 22 24 21 21 21	REC.

COMMENTS:

### WATER SEMIVOLATILE MATRIX SPIKE/MATRIX SPIKE DUPLICATE RECOVERY

SAMPLE CONCENTRATION (ug/L)	CONCENT!		REC #	QC. LIMIT REC. 12-11 27-12 36- 9 41-11 39- 9
	1		·	REC. 12-11 27-12 36- 9 41-11 39- 9
				12-11 27-12 36- 9 41-11 39- 9
				27-12 36- 9 41-11 39- 9
				36- 9 41-11 39- 9
				41-11 39- 9
				39- 9
				23- 9
			/	46-11
				10- 8
				24- 9
				9-10
				26-12
(ug/L)	REC #	RPD #	RPD	REC.
****	33333			12-11
	-			27-12
	<del></del>  -		1	36- 9
	<del></del>  -		1	41-11
	-		1	39- 9
<del></del>	-		42	23- 9
·	-		31	46-11
			50	10- 8
			38	24- 9
	l I _		50	9-10
	MSD CONCENTRATION (ug/L)	CONCENTRATION % REC #	CONCENTRATION % %	CONCENTRATION

	SPIKE	SAMPLE	MS	MS	QC
•	ADDED		CONCENTRATION	•	LIMI
COMPOUND ·	(ug/Kg)	(ug/Kg)	(ug/Kg)	REC #	
*****	*****	=======================================	*******	32223	***
Phenol					26-
2-Chlorophenol					25-1
1,4-Dichlorobenzene					28-:
N-Nitroso-di-n-prop.(1)					41-
1,2,4-Trichlorobenzene_					38-:
4-Chloro-3-methylphenol					26-
Acenaphthene			<u> </u>		31-
4-Nitrophenol		·			11-
			]	l	28-
2,4-Dinitrotoluene					
2,4-Dinitrotoluene Pentachlorophenol					17-
Pentachlorophenol	SPIKE ADDED	MSD CONCENTRATION		QC L	17- 35-
2,4-Dinitrotoluene Pentachlorophenol Pyrene  COMPOUND		1		QC L:	17-: 35-:
PentachlorophenolPyrene	ADDED	CONCENTRATION	8 8	RPD	17- 35- IMIT: REG
PentachlorophenolPyreneCOMPOUND	ADDED	CONCENTRATION	8 8	RPD	IMIT: REC ===: 26- 25-:
COMPOUND Phenol 2-Chlorophenol 1,4-Dichlorobenzene	ADDED	CONCENTRATION	8 8	RPD 35 50 27	IMIT: RE: 26- 25-: 28-:
COMPOUND  Phenol 2-Chlorophenol 1,4-Dichlorobenzene N-Nitroso-di-n-prop.(1)	ADDED	CONCENTRATION	8 8	RPD 35 50 27 38	IMIT: REG 26- 25- 28- 41-
COMPOUND Phenol 2-Chlorophenol 1,4-Dichlorobenzene N-Nitroso-di-n-prop.(1) 1,2,4-Trichlorobenzene	ADDED	CONCENTRATION	8 8	RPD 35 50 27 38 23	IMIT: REC====================================
COMPOUND  Phenol 2-Chlorophenol 1,4-Dichlorobenzene N-Nitroso-di-n-prop.(1) 1,2,4-Trichlorobenzene 4-Chloro-3-methylphenol	ADDED	CONCENTRATION	8 8	RPD 35 50 27 38 23 33	IMIT: REC====================================
COMPOUND  Phenol 2-Chlorophenol 1,4-Dichlorobenzene N-Nitroso-di-n-prop.(1) 1,2,4-Trichlorobenzene 1-Chloro-3-methylphenol Acenaphthene	ADDED	CONCENTRATION	8 8	RPD 35 50 27 38 23 33	IMIT: REC 25- 28- 41- 38- 26- 31-
COMPOUND  Phenol 2-Chlorophenol 1,4-Dichlorobenzene N-Nitroso-di-n-prop.(1) 1,2,4-Trichlorobenzene 1-Chloro-3-methylphenol 1.4-Nitrophenol	ADDED	CONCENTRATION	8 8	RPD 35 50 27 38 23 33 19 50	IMIT: RE: 26- 28- 41- 38- 26- 31- 11-
COMPOUND  Phenol 2-Chlorophenol 1,4-Dichlorobenzene N-Nitroso-di-n-prop.(1) 1,2,4-Trichlorobenzene 4-Chloro-3-methylphenol Acenaphthene 4-Nitrophenol 2,4-Dinitrotoluene	ADDED	CONCENTRATION	8 8	RPD 35 50 27 38 23 33 19 50 47	IMIT: REG 26- 28- 41- 38- 31- 11- 28-
COMPOUND  COMPOUND  Phenol  C-Chlorophenol  1,4-Dichlorobenzene  N-Nitroso-di-n-prop.(1)  1,2,4-Trichlorobenzene  4-Chloro-3-methylphenol  Acenaphthene  4-Nitrophenol	ADDED	CONCENTRATION	8 8	RPD 35 50 27 38 23 33 19 50	IMIT: REC 26- 28- 41- 38- 26- 31- 11-

#### 3E

### WATER PESTICIDE MATRIX SPIKE/MATRIX SPIKE DUPLICATE RECOVERY

	SPIKE	SAMPLE	MS		MS	
COMPOUND	ADDED (ug/L)	CONCENTRATION (ug/L)	CONCENT		REC #	,
	1	, ,			=====	-
gamma-BHC (Lindane)_						_
Heptachlor						
Aldrin						
Dieldrin					<del></del>	
Endrin		<u> </u>				_
			1		l	
4,4'-DDT						_
	SPIKE ADDED	MSD CONCENTRATION		% RPD #	QC L	
4,4'-DDTCOMPOUND	1			i .	_	
COMPOUND gamma-BHC (Lindane)	ADDED (ug/L)	CONCENTRATION	•	RPD #	RPD ====== 15	
COMPOUND  gamma-BHC (Lindane)  Heptachlor	ADDED (ug/L)	CONCENTRATION	•	RPD #	RPD ====== 15 20	
COMPOUND  gamma-BHC (Lindane)  Heptachlor  Aldrin	ADDED (ug/L)	CONCENTRATION	•	RPD #	RPD ====================================	
COMPOUND  gamma-BHC (Lindane)  Heptachlor  Aldrin  Dieldrin	ADDED (ug/L)	CONCENTRATION	•	RPD #	RPD 15 20 22 18	
COMPOUND  gamma-BHC (Lindane)  Heptachlor  Aldrin	ADDED (ug/L)	CONCENTRATION	•	RPD #	RPD ====================================	

#### 3F SOIL PESTICIDE MATRIX SPIKE/MATRIX SPIKE DUPLICATE RECOVERY

Lab	Name:		Contract:	_
Lab	Code:	Case No.:	SAS No.:	SDG No.:
lati	ix Spike - EPA	Sample No.:		

COMPOUND	SPIKE	SAMPLE	MS	MS	QC.
	ADDED	CONCENTRATION	CONCENTRATION	%	LIMITS
	(ug/Kg)	(ug/Kg)	(ug/Kg)	REC #	REC.
gamma-BHC (Lindane) Heptachlor Aldrin Dieldrin Endrin 4,4'-DDT					46-127 35-130 34-132 31-134 42-139 23-134

	SPIKE	MSD	MSD			
	ADDED	CONCENTRATION	*	•	CC L	IMITS
COMPOUND	(ug/Kg)	(ug/Kg)	REC #	RPD #	RPD	REC.
***************	******	*****			======	
gamma-BHC (Lindane)					50	46-127
Heptachlor					31	35-130
Aldrin					43	34-132
Dieldrin					38	31-134
Endrin					45	42-139
4,4'-DDT					50	23-134

# Column to be used to flag recovery and RPD values with an asterisk

*	Va.	lues	outsi	lde	of	QC	lim.	its
---	-----	------	-------	-----	----	----	------	-----

	out ofout	outside li	limits		
COMMENTS:		· · · · · · · · · · · · · · · · · · ·	 	 <u>-</u>	· ——

#### 4A VOLATILE METHOD BLANK SUMMARY

Instrument ID:

	,	VOLATILE MET	OD BLANK SUMMARY	.
Lab	Name:		Contract:	
Lab	Code:	Case No.:	SAS No.:	SDG No.:
Lab	File ID:		Lab Samp	le ID:
Date	Analyzed:		Time Ana	lyzed:

EPA SAMPLE NO.

THIS METHOD BLANK APPLIES TO THE FOLLOWING SAMPLES, MS AND MSD:

	1						
	EPA		LAB		LA		TIME
	SAMPLE		SAMPLE	ID	FILE	ID	ANALYZED
	202###2	*****	********	***	*******	*****	2222######
01				· · · · · · · ·			l
02							
03							
04							
05							
06							
07							
08		-					
09							
10							<del></del>
11							
12							
13		<del></del>					
14	<del></del>		<del></del>		<del></del>		
15							
16							
17							
	<del></del>						
18	<del></del>						
19							
20	<del></del>						
21							
22					<del></del>		
23							<del></del>
24							 
25							
26							
27	•						
28							-
29							
30					<del></del>		
1							

28 29 30			
COMMENTS:	 	 	·
page of _			

# . 4B SEMIVOLATILE METHOD BLANK SUMMARY

ID:			Lab Sample	e ID:
nt I	D:	<del> </del>	Date Extra	acted:
(soi	.l/water)		Date Analy	yzed:
ow/m	ned)	<del>-</del>	Time Analy	yzed:
IS M	ETHOD BLANK A	APPLIES TO THE	FOLLOWING SAMPLI	ES, MS AND MSD:
}	EPA	LAB	LAB	DATE
- 1	SAMPLE NO.	SAMPLE ID	FILE ID	ANALYZED
0.1		*********		********
01				
03				
04				
05				
06				
07		<del></del>	<del></del>	
08	<del></del>			
09	<del></del>			
11	<del></del>	<del></del>		<del></del>
12				
13	<del></del>	······································		
14				
15				
16			 	
17				
18	<del></del>			
20				<del></del>
21		<del></del>	[ <del></del>	<del></del>
22				
23				
24				
25				
26				
27	<del></del>			
28				
29 30				
301			<del></del>	

OLM03.0

# 4C PESTICIDE METHOD BLANK SUMMARY

EPA SAMPLE NO.
----------------

Lab Name:		Contrac	t:		
Lab Code:	Case No.:	SAS No	.:	SDG No.: _	
Lab Sample ID:		_ Lab	File ID:	****	
Matrix:(soil/wat	er)	Extr	action:(Sep	F/Cont/Sonc	)
Sulfur Cleanup:	(Y/N) <u> </u>	Date	Extracted:		<del></del>
Date Analyzed (1	):	Date	Analyzed (	2):	
Time Analyzed (1	):	Time	Analyzed (	2):	<del></del>
Instrument ID (1)	):	Inst	rument ID (	2):	
GC Column (1):	ID:	(mm) GC Co	olumn (2):		ID:(mm)
THIS ME	THOD BLANK AP	PLIES TO THE FO	LLOWING SAM	PLES, MS AN	D MSD:
					!
	EPA	LAB	DATE	DATE	
•	SAMPLE NO.		ANALYZED 1		
01				335555555	
01 02					
03	<u> </u>			<del></del>	
04		<del></del>	<del></del> -	ļ - <del></del>	
05				ļ <del></del>	
06					
07					
08					
. 09					
10				-	
11					
12					
13					
14				ļ	
15					
16					
17					
18   19					•
20					
21		<del></del>		. —————	
22					
23					i
24					
25					
26					
COMMENTS:					
page of				•	

5A

#### VOLATILE ORGANIC INSTRUMENT PERFORMANCE CHECK BROMOFLUOROBENZENE (BFB)

Lab Nar	me:	Contract:	<del></del>	
Lab Cod	de: Case No.:	SAS No.:	SDG No.:	***
Lab Fi	le ID:	BFB Inje	ection Date:_	
Instru	ment ID:	BFB Inje	ection Time:	<del></del>
GC Colu	ımn: ID:(mm)	Heated I	Purge: (Y/N) _	
	ION ABUNDANCE CRITERIA	<del> </del>	Ĭ	RELATIVE ABUNDANCE
	8.0 - 40.0% of mass 95 30.0 - 66.0% of mass 95 Base peak, 100% relative abu			
96	5.0 - 9.0% of mass 95			
174	50.0 - 120.0% of mass 95			( ):
175	4.0 - 9.0 % of mass 174			( ):
	93.0 - 101.0% of mass 174			( ):
177	5.0 - 9.0% of mass 176		1	( )2

1-Value is % mass 174

2-Value is % mass 176

THIS CHECK APPLIES TO THE FOLLOWING SAMPLES, MS, MSD, BLANKS, AND STANDARDS:

	EPA	LAB	LAB	DATE	TIME
	SAMPLE NO.	SAMPLE ID	FILE ID	ANALYZED	ANALYZED
	*********	=======================================		****	
01					
02	<del></del>				
03		<u> </u>			
04					
05					<u></u>
06					
07					
08			<del></del>		
09					
10					
11					
12					
13	<u> </u>				
14					
15					
16					
17					
18					
19					
20					
21					
22					

page \_\_ of \_\_

# SEMIVOLATILE ORGANIC INSTRUMENT PERFORMANCE CHECK DECAFLUOROTRIPHENYLPHOSPHINE (DFTPP)

Lab Nar	me:Contract:	•
Lab Co	de: Case No.: SAS No.: SDG	No.:
Lab Fil	le ID: DFTPP Injection Date	e:
Instru	ment ID: DFTPP Injection Time	e:
m/e	ION ABUNDANCE CRITERIA	RELATIVE ABUNDANCE
51		
68	30.0 - 80.0% of mass 198	]( )1
70	Less than 2.0% of mass 69 25.0 - 75.0% of mass 198	( )1
197	Less than 1.0% of mass 198	
198 199	Base Peak, 100% relative abundance 5.0 to 9.0% of mass 198	
275	10.0 - 30.0% of mass 198	
365	Greater than 0.75% of mass 198	
441	Present, but less than mass 443	
, ,	40.0 - 110.0% of mass 198	
443	15.0 - 24.0% of mass 442	( )2

THIS CHECK APPLIES TO THE FOLLOWING SAMPLES, MS, MSD, BLANKS, AND STANDARDS:

1-Value is % mass 69 2-Value is % mass 442

- 1	EPA	LAB	LAB	DATE	TIME
ŀ	SAMPLE NO.	SAMPLE ID	FILE ID	ANALYZED	ANALYZED
	*****	=======================================		*****	
01		•			
02	<del></del>				
03	<del></del>				
04					
05					
06					
07		<del></del>			
08					<del></del>
09					
10					
11					
12					
13					
14					
15					
16					
17					
18					
19			<del></del>		<del></del>
20 -	<del></del>				<del></del>
21					
22 🛛					

FORM V SV

page \_\_ of \_\_

### VOLATILE ORGANICS INITIAL CALIBRATION DATA

Lab Name:	<del></del>	Contra	ct:				
Lab Code: Case No.:		SAS N	o.:	s	DG No.:		_
Instrument ID: Ca	librati	on Date	(8):				_
Heated Purge: (Y/N) Ca	librati	on Time	s:			<u>-</u>	<del>-</del>
GC Column: ID:	(mm)						
LAB FILE ID: RRF10	E		RRF2	0 =			1
LAB FILE ID: RRF10 RRF50 = RRF100	=		RRF2	00=			
	<del></del>		, <del></del>	<del>,</del>	,	<del>,</del>	l
COMPOUND				RRF100			
	1	)	1	1		1	
Chloromethane	⊦		<b></b>				ļ;
Bromomethane Vinyl Chloride					ļ		<b> </b> ,
Methylene Chloride			<u> </u>				
Acetone							
Carbon Digulfide	ł						
1.1-Dichloroethene	†						;
1,1-Dichloroethene 1,1-Dichloroethane	*					<b></b>	<u>-</u>
1,2-Dichloroethene (total)_							
Chloroform	*				<del></del>	<del></del>	[ <del></del> ;
1,2-Dichloroethane	*						ļ,
2-Butanone							
1,1,1-Trichloroethane	*						,;
Carbon Tetrachloride	*						
Bromodichloromethane	*						
1 2-Dichlerenses	1						[ <del></del>
cis-1,3-Dichloropropene Trichloroethene	*						
Trichloroethene	*						
Dibromochloromethane	*						,
1,1,2-Trichloroethane	*						*
Benzene	*						
trans-1,3-Dichloropropene	*						
Bromoform	*						·
4-Methyl-2-Pentanone							[—— <sub>]</sub>
2-Hexanone							
Tetrachloroethene	*						*
1,1,2,2-Tetrachloroethane	*						
Toluene .	*						*
Chlorobenzene	*			1			*
Ethylbenzene	*						<u>+</u>
Styrene	*			<del>  </del>			*
Xylene (total)	*						
	*=====	=====	=====	=====	*****		=====
Toluene-d8			 				
Bromofluorobenzene	*						*
1,2-Dichloroethane-d4							
			<del></del>	!			

<sup>\*</sup> Compounds with required minimum RRF and maximum %RSD values.
All other compounds must meet a minimum RRF of 0.010.

### SEMIVOLATILE ORGANICS INITIAL CALIBRATION DATA

Lab Name:	<del></del>	Contra	ct:				
Lab Code: Case No.:		SAS N	o.:	s	DG No.:		<del></del>
Instrument ID: Ca	librati	on Date	(8):				<del>_</del>
Ca	librati	on Time	s:	<del> </del>	_		_
			•				
LAB FILE ID: RRF20	<b></b>		RRF50	=		·····	1
RRF80 = RRF120	=		RRF16		•		İ
							.
			ļ				-
COMPOUND	RRF20	RRF50	RRF80	1	RRF160		RSD
	*****	****	*****	*****	*****	*****	
Phenol	<u> </u>					ļ	
District Officer	<u></u>		<del></del>		<del></del>		·
2-Chlorophenol 1,3-Dichlorobenzene	<u> </u>		<del></del>	ļ			·
1,4-Dichlorobenzene	<u> </u>		<del></del>				·
1,2-Dichlorobenzene					<del></del>		·
2-Methylphenol.	<u> </u>		I——				
2,2'-oxybis(1-Chloropropane)	i	<del></del>		1	<del></del>	<del></del>	<del></del>
	<u> </u>		<del></del>			l	
N-Nitrogo-di-n-propylamine		<del></del>	ļ <del></del>				·
Hexachloroethane						<del></del>	ļ <del></del>
Nitrobenzene '	<del>,</del>		·				
Isophorone	.——						
2-Nitrophenol	,——						
2,4-Dimethylphenol	,						-
bis(2-Chloroethoxy)methane	·						
2,4-Dichlorophenol	,						
1,2,4-Trichlorobenzene							
Naphthalene	•						
4-Chloroaniline							
Hexachlorobutadiene							
4-Chloro-3-methylphenol							
2-Methylnaphthalene							
Hexachlorocyclopentadiene							
2,4,6-Trichlorophenol							[
2,4,5-Trichlorophenol*	·						
2-Chloronaphthalene	·						
2-Nitroaniline							<b> </b>
Dimethylphthalate							<b></b>
Acenaphthylene							
2,6-Dinitrotoluene *	' <del></del>						
3-Nitroaniline	,						
Acenaphthene *	·						
2,4-Dinitrophenol					<del></del>		
Dibenzofuran	,						
2.4-Dinitrotoluene	,						

<sup>\*</sup> Compounds with required minimum RRF and maximum %RSD values.
All other compounds must meet a minimum RRF of 0.010.

#### 6C

#### SEMIVOLATILE ORGANICS INITIAL CALIBRATION DATA

Contract:

Lab Name:

nstrument ID: Ca	librati	on Date	(s):	<del></del>	_		_
Ca	librati	on Time	s:		_	<del></del>	_
LAB FILE ID: RRF20	=	<del> </del>	RRF50				1
RRF80 = RRF120		<del></del>	_	0=		_	
			_ ////			<del></del>	_
COMPOUND	RRF20	RRF50	RRF80	RRF120	RRF160	RRF	R
· · ·				****		*====	==
Diethylphthalate							<u> </u>
4-Chlorophenyl-phenylether_	*						
Fluorene	*						
-Nitroaniline							<b> </b>
4,6-Dinitro-2-methylphenol							
N-Nitrosodiphenylamine (1)							
-Bromophenyl-phenylether	*						
iexachlorobenzene	*						
Pentachlorophenol	*						
Phenanthrene	*						l
Inthracene	*						
Carbazole							
i-n-butylphthalate							
luoranthene	*						$  \bot  $
Pyrene	*						
Sutylbenzylphthalate							
3,3'-Dichlorobenzidine							
Benzo(a)anthracene	*						
Chrysene	*						
ois(2-Ethylhexyl)phthalate							
Di-n-octylphthalate							
Benzo(b)fluoranthene	*						
Benzo(k)fluoranthene	*						
	*						
	*		·				
Dibenz(a,h)anthracene	*						
Benzo(g,h,i)perylene	*						_
			======				·
itrobenzene-d5	1				1		l
-Fluorobiphenyl	*						1
erphenyl-d14	*						
henol-d5	*						
-Fluorophenol_	*						
2,4,6-Tribromophenol							
-Chlorophenol-d4	*	<del></del>			<del></del>		
.,2-Dichlorobenzene-d4	*						
·,							

<sup>\*</sup> Compounds with required minimum RRF and maximum %RSD values.
All other compounds must meet a minimum RRF of 0.010.

#### PESTICIDE INITIAL CALIBRATION OF SINGLE COMPONENT ANALYTES

Lab Name:		Contract:	
Lab Code:	Case No.:	SAS No.:	SDG No.:
Instrument ID:	Level (x	low): low mid	high
GC Column:	(mm)	Date(s) Analyzed:	

	RT O	F STAND	ARDS	MEAN	RT W	INDOW
COMPOUND	LOW	MID	HIGH	RT	FROM	TO
alpha-BHC		****		****		
beta-BHC						
delta-BHC						ļ
gamma-BHC (Lindane)						
Heptachlor						
Aldrin						
Heptachlor epoxide						<del></del>
Endosulfan I						
Dieldrin						
4,4'-DDE						
Endrin						
Endosulfan II						
4.4'-DDD						
Endosulfan sulfate						
4,4'-DDT						
Methoxychlor						
Endrin ketone						
Endrin aldehyde						
alpha-Chlordane ·						
gamma-Chlordane				***		
**************	======	======	=====	*****	=====	EE===:
Tetrachloro-m-xylene	İ					
Decachlorobiphenyl						
	<del></del> -					

<sup>\*</sup> Surrogate retention times are measured from Standard Mix A analyses.

Retention time windows are  $\pm$  0.05 minutes for all compounds that elute before Heptachlor epoxide,  $\pm$ 0.07 minutes for all other compounds, except  $\pm$ 0.10 minutes for Decachlorobiphenyl.

#### 6E

### PESTICIDE INITIAL CALIBRATION OF SINGLE COMPONENT ANALYTES

Lab Name:		_ Contrac	t:	_	
Lab Code: Car	se No.:	SAS No	.:	SDG No.:	
Instrument ID:	Level	(x low): lo	w mid	high	
GC Column:	ID:(m	m) Date(s	) Analyzed:		
				٠.	•
	•				
		CALIBRATI	ON FACTORS		
COMPOUND	LOW	MID	HIGH	MEAN	*RSD
•	1	ſ			*****
alpha-BHCbeta-BHC				<del></del>	
delta-BHC	·····				
gamma-BHC (Lindane)	<del></del>				
Heptachlor		<del></del>		<del></del>	
Aldrin					
Heptachlor epoxide		- <del></del>		<del></del>	
Endosulfan I					
Dieldrin					
4,4'-DDE	•				
Endrin		·———			
Endosulfan II					· · ·
עטער אין אין					
Endosulfan sulfate					
4,4'-DDT					
Methoxychlor		·			
Endrin ketone					
Endrin aldehyde					
alpha-Chlordane					
gamma-Chlordane					
	*******	*******	=========	22222222	252322
Tetrachloro-m-xylene					

Decachlorobiphenyl

<sup>\*</sup> Surrogate calibration factors are measured from Standard Mix A analyses.

### PESTICIDE INITIAL CALIBRATION OF MULTICOMPONENT ANALYTES

rap wame: _		Contract:	<del></del>
Lab Code:	Case No.:	SAS No.:	SDG No.:
Instrument	ID:	Date(s) Analyze	ed:
GC Column:	ID:	_(mm)	

***************************************	AMOUNT			RT W	INDOW	CALIBRATION	
COMPOUND	(ng)	PEAK	RT	FROM	TO	FACTOR	
	****	====			****		
Toxaphene		*1	<u> </u>		<b> </b>		
		*2 *3				<del></del>	
		4					
•		5	<del></del>				
Aroclor 1016	<del></del>						
1200001 2020		*2				· — —	
		*3					
		4					
		5_					
Aroclor 1221		*1					
		*2					
		*3		 			
	•	4				<del></del>	
<del></del>		5-					
Aroclor 1232		*1					
		*2					
		*3					
		<b>4</b> 5					
Aroclor 1242						<del></del>	
RIOCIOI 1242		*2					
		*3					
		4				<del></del>	
		5					
Aroclor 1248		*1					
		*2					
·		*3					
		4					
		_ 5					
Aroclor 1254		*1					
		*2					
	Ì	*3					
ſ	. 1	4					
		5_				<del></del>	
Aroclor 1260		*1					
		*2					
j	ļ	*3					
	ļ	4 5				<del></del>	
		٥				<del></del>	
				l			

<sup>\*</sup> Denotes required peaks

# 6G PESTICIDE ANALYTE RESOLUTION SUMMARY

Lab Name: Con	ontract:	
Lab Code: Case No.: S	AS No.: SDG No.:	
GC Column (1): ID: (mm)	Instrument ID (1):	
EPA Sample No. (Standard 1):	Lab Sample ID (1):	-
Date Analyzed (1):	Time Analyzed (1):	
	RESOLUTION	
ANALYTE	RT (%)	
01 02	_	
03	_	
04		
05		
06		
07		
08		
09		
GC Column (2): ID:(mm)	Instrument ID (2):	
EPA Sample No. (Standard 2):	Lab Sample ID (2):	-
Date Analyzed (2):	Time Analyzed (2):	
	RESOLUTION	
ANALYTE	RT (%)	
01	ま   日本本に正正   日本社会を単独事を書	
02		
03	-	
04	-	
05	<u> </u>	
.06		
07		
.08		

# 6H PERFORMANCE EVALUATION MIXTURE (PEM)

Lab Name:		Contract:_	ntract:			
Lab Code:	Case No.:	SAS No.: _	SDG No.:			
GC Column (1)	: ID:(mm	) Instrume	nt ID (1):			
EPA Sample No. (Standard 1): Lab Sample ID (1):						
Date Analyzed (1): Time Analyzed (1):						
ı	<del></del>		RESOLUTION			
	ANALYTE	RT	(%)			
	****************	1 }				
01		.				
02						
03						
04		.				
05						
06   07		-				
. 08		·	<del></del>			
001	<del></del>	I	l			
GC Column (2): ID: (mm) Instrument ID (2):						
EPA Sample No. (Standard 2): Lab Sample ID (2):						
Date Analyzed (2): Time Analyzed (2):						
1			RESOLUTION			
	ANALYTE	RT	(%)			
01		.				
02						
03		.				
04						
05	<del></del>	-				
06	<del></del>					
0/		.				

#### 6I INDIVIDUAL STANDARD MIXTURE A

Lab Name:	Contract:			
Lab Code:	Case No.:	SAS No.:	SDG No.:	
	ID:(mm		•	
Date Analyzed (	1):	Time Anal	yzed (1):	
	ANALYTE	RT	ESOLUTION (%)	
01 02 03				
04 05 06 07				
08 09 10				
11	<u> </u>	ll		
	ID:(mm			
	2):		yzed (2):	
	ANALYTE	RT	ESOLUTION (%)	
01 02 03 04				
05 06 07				
08 09 10 11				

# 6J INDIVIDUAL STANDARD MIXTURE B

Lab Name:	Contract:					
Lab Code:	Case No.: SA	AS No.:	SDG No.:			
GC Column (1	):(mm)	Instrum	ent ID (1):			
EPA Sample No. (Standard 1):			Lab Sample ID (1):			
Date Analyzed (1): Time Analyzed (1):						
	1	T	RESOLUTION			
	ANALYTE	RT	(%)			
01			l '			
02		· <del> </del>				
03						
04						
05						
06						
07		.				
08		·				
09 10	•	<b> </b>				
11		ļ				
12		ļ				
13						
		' <del></del>	'			
GC Column (2)	):(mm)	Instrume	ent ID (2):			
EPA Sample No. (Standard 2): Lab Sample ID (2):						
Date Analyzed (2): Time Analyzed (2):						
į		I	RESOLUTION			
	ANALYTE	RT	(%)			
	***************	=====				
01		<u> </u>				
02						
03						
04						
05						
06 07						
08						
09						
10	<del></del>					
11						
12						
13	<del></del>					

### VOLATILE CONTINUING CALIBRATION CHECK

Lab Name:		Contract:		
Lab Code:	Case No.: _	SAS No.:	SDG	No.:
Instrument ID:	Cal	ibration Date:_	Tin	ne:
Lab File ID:	Ini	t. Calib. Date(	s):	
Heated Purge: (Y/N)	Ini	t. Calib. Times	:	
GC Column:	ID:	(mm)		

			MIN		MAX
COMPOUND	RRF	RRF50	RRF	%D	&D
Chloromethane	780 <b>5</b> 8		=====		
Bromomethane			0.100		25.0
Vinyl Chloride			0.100		25.0
Chloroethane			0.100		23.0
Methylene Chloride					
Acetone			1		İ
Carbon Disulfide					
1,1-Dichloroethene			0.100		25.0
1,1-Dichloroethane	<del></del>		0.200		25.0
1,2-Dichloroethene (total)			0.200		123.0
Chloroform			0.200		25.0
1,2-Dichloroethane			0.100		25.0
2-Butanone		<del></del>	0.100	<del></del>	23.0
1,1,1-Trichloroethane			0.100		25.0
Carbon Tetrachloride			0.100		25.0
Bromodichloromethane			0.200		25.0
1,2-Dichloropropane			0.200		23.0
cis-1,3-Dichloropropene			0.200		25.0
Trichloroethene			0.300	<del></del>	25.0
Dibromochloromethane			0.100		25.0
1,1,2-Trichloroethane			0.100		25.0
Benzene			0.500		25.0
trans-1,3-Dichloropropene			0.100		25.0
Bromoform			0.100		25.0
4-Methyl-2-Pentanone					
2-Hexanone					
Tetrachloroethene			0.200		25.0
1,1,2,2-Tetrachloroethane	·		0.300		25.0
Toluene			0.400		25.0
Chlorobenzene			0.500		25.0
Ethylbenzene			0.100		25.0
Styrene			0.300		25.0
Xylene (total)			0.300		25.0
- , , , , , , , , , , , , , , , , , , ,	=====	*****	=====	RESEER	====
Toluene-d8					
Bromofluorobenzene			0.200		25.0
1,2-Dichloroethane-d4					

All other compounds must meet a minimum RRF of 0.010.

### SEMIVOLATILE CONTINUING CALIBRATION CHECK

Lab Name:	Contract:	
Lab Code: Case	No.: SAS No.: SDG No.:	
Instrument ID:	Calibration Date: Time:	
Lab File ID:	Init. Calib. Date(s):	
	Init. Calib. Times:	

			MIN		MA
COMPOUND	RRF	RRF50	RRF	₹D	*D
Phenol			0.800		25.0
bis(2-Chloroethyl)ether			0.700		25.0
2-Chlorophenol			0.800		25.0
1,3-Dichlorobenzene			0.600		25.0
1,4-Dichlorobenzene			0.500		25.
1,2-Dichlorobenzene_			0.400		25.
2-Methylphenol			0.700		25.
2,2'-oxybis(1-Chloropropane)					•
4-Methylphenol			0.600		25.
N-Nitroso-di-n-propylamine			0.500		25.0
Hexachloroethane			0.300		25.0
Nitrobenzene			0.200		25.0
Isophorone			0.400		25.0
2-Nitrophenol			0.100		25.0
2,4-Dimethylphenol			0.200		25.0
bis(2-Chloroethoxy)methane			0.300		25.0
2,4-Dichlorophenol			0.200		25.0
1,2,4-Trichlorobenzene			0.200		25.0
Naphthalene			0.700		25.0
4-Chloroaniline					
Hexachlorobutadiene					
4-Chloro-3-methylphenol			0.200		25.0
2-Methylnaphthalene			0.400		25.0
Hexachlorocyclopentadiene					
2,4,6-Trichlorophenol			0.200		25.0
2,4,5-Trichlorophenol			0.200		25.0
2-Chloronaphthalene			0.800		25.0
2-Nitroaniline					
>imethylphthalate					
Acenaphthylene			0.900		25.0
2,6-Dinitrotoluene			0.200		25.0
3-Nitroaniline					
Acenaphthene			0.900		25.0
2,4-Dinitrophenol					
-Nitrophenol			1		
Dibenzofuran			0.800		25.0
2,4-Dinitrotoluene			0.200		25.0
-,					

All other compounds must meet a minimum RRF of 0.010.

#### 7C

#### SEMIVOLATILE CONTINUING CALIBRATION CHECK

Lab Name:	Contract:	
Lab Code: Cas	se No.: SAS No.:	SDG No.:
Instrument ID:	Calibration Date:	Time:
Lab File ID:	Init. Calib. Date(s):	
	Init. Calib. Times:	

	l		MIN		MAX
COMPOUND	RRF	RRF50	RRF	%D	<b>%</b> D
	=====	=====	2222		
Diethylphthalate					
4-Chlorophenyl-phenylether			0.400		25.0
Fluorene	<b> </b>	<b></b>	0.900		25.0
4-Nitroaniline					
4,6-Dinitro-2-methylphenol	ļ				
N-Nitrosodiphenylamine (1)					
4-Bromophenyl-phenylether	<b> </b>		0.100		25.0
Hexachlorobenzene			0.100		25.0
Pentachlorophenol	<b> </b>		0.050		25.0
Phenanthrene			0.700		25.0
Anthracene	.   <del></del>		0.700		25.0
Carbazole					1
Di-n-butylphthalate					Ì
Fluoranthene			0.600		25.0
Pyrene			0.600		25.0
Butylbenzylphthalate					1
3,3'-Dichlorobenzidine					
Benzo(a)anthracene			0.800		25.0
Chrysene			0.700		25.0
bis(2-Ethylhexyl)phthalate			-		
Di-n-octylphthalate					]
Benzo(b) fluoranthene			0.700		25.0
Benzo(k)fluoranthene			0.700		25.0
Benzo(a)pyrene			0.700		25.0
Indeno(1,2,3-cd)pyrene			0.500		25.0
Dibenz(a,h)anthracene			0.400		25.0
Benzo(g,h,i)perylene			0.500		25.0
			****	=====	8822
Nitrobenzene-d5		'	0.200		25.0
2-Fluorobiphenyl		-	0.700		25.0
Terphenyl-d14			0.500		25.0
Phenol-d5			0.800		25.0
2-Fluorophenol			0.600		25.0
2,4,6-Tribromophenol					<del>-</del>
2-Chlorophenol-d4			0.800		25.0
1,2-Dichlorobenzene-d4	<del></del>		0.400		25.0

<sup>(1)</sup> Cannot be separated from Diphenylamine

All other compounds must meet a minimum RRF of 0.010.

# 7D PESTICIDE CALIBRATION VERIFICATION SUMMARY

Lab Name:		Conti	ract:			
Lab Code: Case No.:		SAS	No.: _	sdc	No.:	<del></del>
GC Column: ID:	(mm)	Init.	Calib.	Date(s):		
EPA Sample No.(PIBLK):			Dat	e Analyzed	l :	
Lab Sample ID (PIBLK):			Tim	e Analyzed	l :	
EPA Sample No. (PEM):	<del></del>		Dat	e Analyzed	ı :	
Lab Sample ID (PEM):			Tim	e Analyzed	l :	
PEM		RT V	MODUIN	CALC	NOM	
COMPOUND	RT	FROM	TO	AMOUNT (ng)	AMOUNT (ng)	<b>%</b> D
			·	_	<b>尼亚亚联巴拉拉克</b>	
alpha-BHC beta-BHC						
gamma-BHC (Lindane)			-			<del></del>
Endrin						
4,4'-DDT						
Methoxychlor						
			.			
4,4'-DDT % breakdown (1):		E	indrin %	breakdown	(1):	<del></del>

# 7E PESTICIDE CALIBRATION VERIFICATION SUMMARY

Lab Name:		contr	act:			
Lab Code: Case No	o.:	SAS	No.:	SD0	3 No.:	
GC Column: ID:	( mm )	Init.	Calib.	Date(s):_	<del></del>	
EPA Sample No.(PIBLK):			Date	e Analyze	d :	<del></del>
Lab Sample ID (PIBLK):	<u>-</u>		Time	e Analyze	d :	
EPA Sample No.(INDA):			Date	e Analyze	d :	
Lab Sample ID (INDA):			Time	e Analyze	d :	<del></del>
INDIVIDUAL MIX A		RT W	INDOW	CALC	NOM	
COMPOUND	RT	FROM		AMOUNT	AMOUNT	8D
		}	1	(ng)	(ng)	
	=======	======	*****			3==3#
alpha-BHC				}	ļ	
gamma-BHC (Lindane)		1				
Heptachlor	-					
Heptachlor Endosulfan I	-	<del></del>		<del></del>		
Dieldrin	-					
Endrin	_	<del></del>	i i			
		<u> </u>				
A A'-DDT	1	1				
Methoxychlor	•	1	1			
Tetrachloro-m-xylene	-					
Decachlorobiphenyl	_					
EPA Sample No.(INDB): Lab Sample ID (INDB):					i : i :	
INDIVIDUAL MIX B		RT W	INDOW	CALC	МОМ	
COMPOUND	RT	FROM	TO	AMOUNT	THUOMA	%D
		]	]	(ng)	(ng)	
		######	=====			=====
beta-BHC						
delta-BHC				<del></del>		
Aldrin	_			<del></del>		
Heptachlor epoxide	_					
4,4'-DDE	-					
Endosulfan II	-					
Endosulfan sulfate					i <del></del>	
Endrin ketone					, <del></del>	
Endrin aldehyde	_					
alpha-Chlordane	-		<del></del>	<del></del>	, <del></del>	
gamma-Chlordane				<del></del>		<del></del>
Tetrachloro-m-xylene				<del></del>		
Decachlorobiphenyl						
Pacacutotopthuanit -						
			<b></b>			

### 8A VOLATILE INTERNAL STANDARD AREA AND RT SUMMARY

Lab	Code:	Case No	.:	SAS No.:	<del></del>	SDG No.:	<del></del>
Lab	File ID (Sta	ndard):	·		Date A	nalyzed:	
Inst	trument ID:				Time A	nalyzed:	
GC (	Column:	ID:	(mm)		Heated	Purge: (Y/1	4)
ĺ		IS1(BCM)		IS2(DFB)		IS3(CBZ)	
		AREA #	}				
	12 HOUR STD			*****	=======	*******	
	UPPER LIMIT			<del></del>	<u> </u>	<del></del>	
	LOWER LIMIT						<del></del>
	****	******	======		======	*********	*=====
	EPA SAMPLE NO.						
	********	*****		*********	*****	********	*=====
01						<del></del>	
02						<del></del>	
04							
05							
06				· · ·			
07							
80							
09							
11							
12				<del></del>			
13							
14							
15							
16							
18		<del></del>					
19							
20							
21							
22		}					

OLM03.0

page \_\_ of \_\_

#### SEMIVOLATILE INTERNAL STANDARD AREA AND RT SUMMARY

Lab Name: Contract:

				5A5 NO.: _		SDG No.:	
	le ID (Stand	ard):		1	Date Ana	lyzed:	
stru	ment ID:			;	Time Ana	lyzed:	
1-	·	IS1(DCB)	T	IS2(NPT)		IS3(ANT)	
	•	AREA #	,		RT #		RT
-			,		1	]	1
.	12 HOUR STD	1					}
- 1	UPPER LIMIT				·		
1	LOWER LIMIT						
-	*========	######################################	222222	**======	****	========	*****
	EPA SAMPLE NO.			-			
}-	******	******	*=====	*****	======	353853E3E2	=====
01		1	1			!	
02					•		
03	<del></del>						
04							
05						·	
06							
07							
08]				•			
09 _							
10							
11							
12							
13			l				
14 _			l				
15 _			ļ				
16			l				
17 _							
18 _			<b></b>				
19 _			· 				
20							
21							
22 _			i				

#### 8c '

#### SEMIVOLATILE INTERNAL STANDARD AREA AND RT SUMMARY

Lab N	ame:			Contract:		_	
Lab C	ode:	Case No.:		SAS No.: _		SDG No.:	
Lab F	ile ID (Stand	ard):		1	Date Ana	lyzed:	<del></del>
Instr	ument ID:	<del></del>			Time Ana	lyzed:	
	[	IS4(PHN)		IS5(CRY)		IS6(PRY)	
		AREA #	RT #	AREA #	RT #	AREA #	RT #
	*********			*****	*****	******	
	12 HOUR STD						
	UPPER LIMIT		·				
	LOWER LIMIT	}					ļ
			222222		======		
	EPA SAMPLE			,			
	NO.				,	•	
	*********	=======================================			*****	*****	22222
01							
02							
03	]		]				
04							
05							
06	<del></del>						
07	<del></del>						<del></del>
08							
09	<del></del>						
10							
11	<del></del>						<u> </u>
12 13							
14	<del></del>						
15							
16							
17							
18		<del></del>		<del></del>			
19							
20							<del></del>
21							<del></del>
22							<del></del>
IS IS AF AF RT	54 (PHN) = Pho 55 (CRY) = Chr 56 (PRY) = Per REA UPPER LIMI T UPPER LIMIT T LOWER LIMIT Column used to Values outside	rysene-d12 rylene-d12 T = +100% o T = - 50% o = +0.50 min = -0.50 min	of intermof intermotors of intermotors of outes of ernal sta	al standard internal st	area andard R andard R	T	risk.
page _	of ,						

# 8D PESTICIDE ANALYTICAL SEQUENCE

		PESTICIDE	ANALYTICAL	SEQUENCE		
Lab Name:			Contract	:	_	
Lab Code:	: C	ase No.:	SAS No.	:	SDG No.:	
GC Column	n:	ID:(mm	) Init. Cal	ib. Date(s)	:	
Instrumer	nt ID:					
THE ANA	ALYTICAL SEQUE SAMPLES	ENCE OF PERFO			URES, BLAI	NKS,
		GATE RT FROM DCB:		IBRATION		
	EPA	LAB SAMPLE ID	DATE	TIME	TCX	DCB
0.1	EZ#8EZEFE	22500220000000000000000000000000000000	*******		******	2322233
01 02				<del></del>	<del></del>	
03						<del></del>
04						
05	<u> </u>					
06						
07						
80 90		<del></del>			ļ <del></del>	ļ ————
10	<u></u>			<del></del>	<del></del>	<del></del>
11	<del></del>					
12						
13						
14		<del></del>		<u> </u>		<del></del>
15 16	<del></del>	<del></del>	<del></del>	ļ <del></del>	ļ	<del></del>
17				J <del></del>	<del></del>	
18						
19						
20						
21						
22						
23 24						
25	<del></del>	<del></del>		<del></del>	<del></del>	
26						
27						
28						
29						
30			<del></del>			
31 32					<del></del>	
22				l		

QC LIMITS

TCX = Tetrachloro-m-xylene ( $\pm$  0.05 MINUTES)

DCB = Decachlorobiphenyl  $(\pm 0.10 \text{ MINUTES})$ 

# Column used to flag retention time values with an asterisk.

\* Values outside of QC limits.

page \_ of \_

### 9A . PESTICIDE FLORISIL CARTRIDGE CHECK

Lab Name:	Co	ontract:	<b></b>		
Lab Code:	Case No.:	SAS No.:	SDG No.:		
Florisil Cartridge	Lot Number:	Date of Analys	is:	<del></del>	
GC Column(1):	ID:(mm)	GC Column(2):		ID:	_(mm)

SPIKE	SPIKE		
ADDED	RECOVERED	8	δc
(ng)	(ng)	REC #	LIMITS
========	=======	=====	2522222
			80-120
			80-120
			80-120
			80-120
			80-120
	•		80-120
			80-120
			80-120
			80-120
			80-120
			80-120
	ADDED (ng)	ADDED RECOVERED (ng)	ADDED RECOVERED % (ng) (ng) REC #

- # Column to be used to flag recovery with an asterisk.
- \* Values outside of QC limits.

THIS CARTRIDGE LOT APPLIES TO THE FOLLOWING SAMPLES, BLANKS, MS, AND MSD:

EPA	LAB	DATE	DATE
SAMPLE NO.	SAMPLE ID	ANALYZED 1	ANALYZED 2
********	202222222	=========	****
01			
02			
03			
04			
05			
06			
07			
08			
09			
10			
11		<del></del>	
12			
13			
14			
15			
16			
17			
18			
19			
20			
21			
22			
23			

page \_\_ of \_\_

#### 9B PESTICIDE GPC CALIBRATION

Lab Name:	Contrac	ct:			
Lab Code: Case No.:	SAS No	o.:	_ SDG 1	No.:	
GPC Column:	Calib	ration Date	e:		
GC Column(1): ID:	(mm) GC (	Column(2):		ID:	(mm)
COMPOUND	SPIKE ADDED (ng)	SPIKE RECOVERED (ng)	REC #	QC. LIMITS REC.	
gamma-BHC (Lindane) Heptachlor Aldrin Dieldrin Endrin 4,4'-DDT	**********			80-110 80-110 80-110 80-110 80-110 80-110	

THIS GPC CALIBRATION APPLIES TO THE FOLLOWING SAMPLES, BLANKS, MS AND MSD:

	EPA	LAB	DATE_	DATE
	SAMPLE NO.	SAMPLE ID	ANALYZED 1	ANALYZED 2
	224252222	7222222222222		*******
01				
02				
03				
04				
05		·		<u></u>
06				
07				
80				
09				
10				
11				
12				<del></del>
13				
14			<del></del> -	
15	<del></del>		<del></del>	
16			<del></del>	
17				
18				
19	<del></del>	<del></del>		
20				
21				
22				
		- <del></del>		
23				
24				
25				
26			<u></u>	

page \_\_ of \_\_

<sup>#</sup> Column to be used to flag recovery values with an asterisk
\* Values outside of QC limits

10A

EPA SAMPLE NO.

#### PESTICIDE IDENTIFICATION SUMMARY FOR SINGLE COMPONENT ANALYTES

	FOR	JINC	JDD COM	-ONENI P	TINDI I DE	<b>'</b>	
ab Name:			Co	ntract:			
ab Code:	_ Case No.	:	s	AS No.:	<del></del>	_ SDG No.:	
ab Sample ID : _				Date(s	) Analy	zed:	
nstrument ID (1)	•			Instru	ment ID	(2):	
Column(1):	ID:		(mm)	GC Col	umn(2):	ID	:
ANALYTE		COL	7,7	I .	INDOW	CONCENTRATION	
	232222222		RT	FROM	TO	CONCENTRATION	i -
		1					
		2		·			
·		1					
		2					
<b>V</b>		1					
		2					
·		1					
		2					-
		1					
		2					
		1					
		2					
	•	1					
		2					
		1					
<del></del>							

page \_ of \_

#### 10B

EPA SAMPLE NO.

# PESTICIDE IDENTIFICATION SUMMARY FOR MULTICOMPONENT ANALYTES

Lab Name:	Co	ntract:	_	
Lab Code:	Case No.: S	AS No.:	SDG No.:	
Lab Sample ID :		Date(s) Analyzed	i:	
Instrument ID (1):		Instrument ID (2	2):	
GC Column(1):	ID: (mm)	GC Column(2):	ID:	_ ( mm )

ANALYTE PEAK RT RT WINDOW CONCENTRATION &D  1 2 3 3 4 5 5 4 5 5 6 6 6 6 6 6 6 6 6 6 6 6 6 6							•	
COLUMN 1		T		RT W	INDOW		MEAN	
COLUMN 1	ANALYTE	PEAK	RT	FROM	TO	CONCENTRATION	CONCENTRATION	%D
COLUMN 1	=======================================	2227	=====	=====	=====	=======================================	==========	*=====
COLUMN 1		1						
COLUMN 1	1:	1						
COLUMN 2 4 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5								
COLUMN 2 4 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5	COLUMN 1	4						
COLUMN 2		5						l
COLUMN 2		1			<del></del>			
COLUMN 2		1			ĺ			
COLUMN 2								
COLUMN 2								
COLUMN 1	COLUMN 2							
COLUMN 1	0020:2: Z	-						
COLUMN 1		Į –	======	35555	=====	*********		
COLUMN 1		j				,		]
COLUMN 1		1				<del></del>		
COLUMN 1						<del></del>		
COLUMN 2	COLUMN 1							
COLUMN 2	COLUMN							
COLUMN 2		3						
COLUMN 2	· :							
COLUMN 2				<del></del>				
COLUMN 2 4 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5								
COLUMN 1		1						
COLUMN 1	COLUMN 2	i -			<del></del>			
COLUMN 1		j -						
COLUMN 1 2 3	=======================================		22222	=====	=====	====================================		*25555
COLUMN 1								
COLUMN 1								
5		ł						
COLUMN 2 4	COLUMN 1	i -						
COLUMN 2 4		5						
COLUMN 2 4								
COLUMN 2 3 4		1						
COLUMN 2 4								
(		3						
5	COLUMN 2	4						

At least 3 peaks for each column are required for identification of multicomponent analytes

page \_\_ of \_\_

#### SAMPLE LOG-IN SHEET

Lab	Name					Page of
Rec	eived by (Print Name)				<del></del>	Log-in Date
Rec	eived By (Signature)	·- · · · · · · · · · · · · · · · · · ·		<del></del>		**************************************
Cas	Number		Sample Delive	ery Group No.		SAS Number
R ema	arks:			Corre	sponding	Remarks:
		•	EPA Sample #	Sample Tag #	Assigned Lab #	Condition of Sample Shipment, etc.
ı	Custody Seal(s)	Present/Absent* Intact/Broken				
2.	Custody Seal Nos.					
3.	Chain-of Custody Records	Present/Absent*				
4.	Traffic Reports or Packing Lists	Present/Absent*				
5.	Airbill	Airbill/Sticker Present/Absent*		·,		
6.	Airbill No.					
7.	Sample Tags	Present/Absent*				
	Sample Tag Numbers 1	Listed/Not Listed on Chain-of-Custody				
1.	Sample Condition	Intact/Broken*/Leaking				
9.	Does information on custody records, traffic reports, and sample tags agree?	Yes/No*				
10.	Date Received at Lab	<del></del>				
11.	Time Received					
	Sample Transfer	·				
	tion ·	Fraction				
Area		Area #	•			
Ву		Ву				
מכ		Ch .				
	Contact \$40 and attach record of r	esolution.				
	ewed By			Logbook No.		
Date				Logbook Page A	10.	

FORM DC-1

OLM03.0

CITY/STATE						
CASE NO.	SDG NO. SDG SDG	NOS. TO	FOLL	wc		
	ivered in the Complete SDG File		orig	inal do	cuments	
		,	DACE	NOs		CHECK
			ROM	-	LAB	EP
nventory Sheet (Fo	orm DC-2) (Do not number)	_				
DG Case Narrative	, ,	_				
DG Cover Sheet/Tra	affic Report	_				
olatiles Data					<del></del>	
. QC Summary						
System Monitori	ng Compound Summary (Form II VO	A)				
Matrix Spike/Ma	trix Spike Duplicate Summary					
(Form III	VOA)				·	
Method Blank Su	mmary (Form IV VOA)					
-	t Performance Check (Form V VOA)	)				
Internal Standa (Form VIII	rd Area and RT Summary VOA)	· _			<del></del>	<del></del>
. Sample Data		_				
TCL Results - (	Form I VOA)					
Tentatively Ide	ntified Compounds (Form I VOA-T	IC)				
Reconstructed to each sample	otal ion chromatograms (RIC) for	c				
For each sample	•					
	nd background-subtracted mass target compounds identified					-
Quantitation	reports `					
Mass spectra best libra	of all reported TICs with three ry matches					
. Standards Data	(All Instruments)	_			•	
Initial Calibra	tion Data (Form VI VOA)					
	eports for all Standards					
	bration Data (Form VII VOA)					
_	tation Reports for all Standards	3				
. Raw QC Data						
BFB	·	_				
Blank Data		_				
Matrix Spike/Ma	trix Spike Duplicate Data					

CASE NO. SDG NO. SDG NOS. TO FOLLOW

Ŀ	SAS NO.	<del></del>	
		PAGE NOS	CHECK
		FROM TO	LAB EP.
5. <u>S</u> e	emivolatiles Data		
а.	QC Summary		
	Surrogate Percent Recovery Summary (Form II SV)		*******
	MS/MSD Summary (Form III SV)	<del></del>	
	Method Blank Summary (Form IV SV)		
	GC/MS Instrument Performance Check (Form V SV)		·
	Internal Standard Area and RT Summary		
	(Form VIII SV)	<del></del>	<del></del>
b.	Sample Data		
	TCL Results (Form I SV-1, SV-2)		
	Tentatively Identified Compounds (Form I SV-TIC)		<del></del>
	Reconstructed total ion chromatograms (RIC) for		
	each sample	•,	
	For each sample:		
	Raw spectra and background-subtracted mass		•
	spectra of target compounds		
	Quantitation reports		<del></del>
	Mass spectra of TICs with three best library matches		
	GPC chromatograms (if GPC performed)		
c.	Standards Data (All Instruments)		
	Initial Calibration Data (Form VI SV-1, SV-2)	•	
	RICs and Quan Reports for all Standards		
	Continuing Calibration Data (Form VII SV-1, SV-2)		
	RICs and Quantitation Reports for all Standards		
	Semivolatile GPC Calibration Data-UV detector traces		
d.	Raw QC Data		
	DFTPP	<del></del>	<del></del>
	Blank Data		<del></del>
	Matrix Spike/Matrix Spike Duplicate Data		
e.	Raw GPC Data	<del></del>	
. Pe	sticides		
a.	QC Summary		
	Surrogate Percent Recovery Summary (Form II PEST)	<del></del>	<del></del>
	MS/MSD Duplicate Summary (Form III PEST)		
	Method Blank Summary (Form IV PEST)		
	•		

FORM DC-2-2 OLM03.0

CASE NO SDG NO SAS NO	SDG NOS. T	o Follo	w	<del></del>	
				<u> </u>	
		PAGE	NOs		CHECK
		FROM	TO	LAB	EP
,				•	
esticides (Cont.)					
. Sample Data			·		
TCL Results - Organic Analysis Data Sheet (Form I PEST)					
Chromatograms (Primary Column)					
Chromatograms from second GC column confirm	mation				
GC Integration report or data system printe					
Manual work sheets			•		
For pesticides/Aroclors confirmed by GC/MS	,				
copies of raw spectra and copies of back subtracted mass spectra of target compou	ground-				
(samples & standards)					
. Standards Data					
Initial Calibration of Single Component Ana (Form VI PEST-1 and PEST-2)	alytes				
Initial Calibration of Multicomponent Analy (Form VI PEST-3)	ytes				<u> </u>
Analyte Resolution Summary (Form VI PEST-4)	)				
Performance Evaluation Mixture (Form VI PES	ST-5)			·	
Individual Standard Mixture A (Form VI PEST	r-6)				
Individual Standard Mixture B (Form VI PEST	T-7)				
Calibration Verification Summary (Form VII PEST-1)					
Calibration Verification Summary (Form VII PEST-2)					
Analytical Sequence (Form VIII PEST)			•		
Florisil Cartridge Check (Form IX PEST-1)					
Pesticide GPC Calibration (Form IX PEST-2)					
Pesticide Identification Summary for Single Component Analytes (Form X PEST-1)	9 ,				
Pesticide Identification Summary for Multicomponent Analytes (Form X PEST-2)				•	
Chromatograms and data system printouts					
A printout of retention times and correspeak areas or peak heights	ponding				
Pesticide GPC calibration data - UV detector traces	or				
Part of Parts					
. Raw QC Data Blank Data					
Matrix Spike/Matrix Spike Duplicate Data			<del> </del>		

PAGE NOS FROM TO LAB  6. Pesticides (Cont.) e. Raw GPC Data  f. Raw Florisil Data  7. Miscallaneous Data Original preparation and analysis forms or copies of preparation and analysis logbook pages Internal sample and sample extract transfer chain- of-custody records Screening records All instrument output, including strip charts from screening activities (describe or list)  6. EPA Shipping/Receiving Documents Airrills (No. of shipments _) Chain-of-Custody Records Sample Tags Sample Tags Sample Tags Sample Log-In Sheet (Lab & DC1) Miscellaneous Shipping/Receiving Records (describe or list)  Internal Lab Sample Transfer Records and Tracking Sheets (describe or list)  Telephone Communication Log	- }	CASE NO.	SDG NO.	SDG NOS.	TO FOLLOW	
6. Pesticides (Cont.) e. Raw GPC Data  f. Raw Florisil Data  7. Miscellaneous Data Original preparation and analysis forms or copies of preparation and analysis logbook pages Internal sample and sample extract transfer chain- of-custody records Screening records All instrument output, including strip charts from screening activities (describe or list)  DPA Shipping/Receiving Documents Airtills (No. of shipments) Chain-of-Custody Records Sample Log-In Sheet (Lab & DC1) Miscellaneous Shipping/Receiving Records (describe or list)  Internal Lab Sample Transfer Records and Tracking Sheets (describe or list)	_ا		3A3 NO			
FROM TO LAB  6. Pesticides (Cont.) e. Raw GPC Data  f. Raw Florisil Data  7. Miscellaneous Data Original preparation and analysis forms or copies of preparation and analysis logbook pages Internal sample and sample extract transfer chain- of-custody records Screening records All instrument output, including strip charts from screening activities (describe or list)					DACE NOS	CHECK
e. Raw GPC Data  f. Raw Florisil Data  / Miscellaneous Data Original preparation and analysis forms or copies of preparation and analysis logbook pages Internal sample and sample extract transfer chain- of-custody records Screening records All instrument output, including strip charts from screening activities (describe or list)		٠				
e. Raw GPC Data  f. Raw Florisil Data  Miscellaneous Data Original preparation and analysis forms or copies of preparation and analysis logbook pages Internal sample and sample extract transfer chain- of-custody records Screening records All instrument output, including strip charts from screening activities (describe or list)  EPA Shipping/Receiving Documents Airtills (No. of shipments) Chain-of-Custody Records Sample Tags Sample Log-In Sheet (Lab & DC1) Miscellaneous Shipping/Receiving Records (describe or list)  Internal Lab Sample Transfer Records and Tracking Sheets (describe or list)	5. P	'esticides (Cont	.)			
Miscellaneous Data Original preparation and analysis forms or copies of preparation and analysis logbook pages Internal sample and sample extract transfer chain- of-custody records Screening records All instrument output, including strip charts from screening activities (describe or list)	_				<del></del>	
Original preparation and analysis forms or copies of preparation and analysis logbook pages  Internal sample and sample extract transfer chain- of-custody records  Screening records  All instrument output, including strip charts from screening activities (describe or list)	f	. Raw Plorisil	Data		· ·	
preparation and analysis logbook pages Internal sample and sample extract transfer chain- of-custody records Screening records All instrument output, including strip charts from screening activities (describe or list)						•
Internal sample and sample extract transfer chain- of-custody records  Screening records All instrument output, including strip charts from screening activities (describe or list)	0:					
Screening records All instrument output, including strip charts from screening activities (describe or list)	I	nternal sample a	and sample extract tran	•	<del></del>	
All instrument output, including strip charts from screening activities (describe or list)  EPA Shipping/Receiving Documents Airbills (No. of shipments) Chain-of-Custody Records Sample Tags Sample Log-In Sheet (Lab & DC1) Miscellaneous Shipping/Receiving Records (describe or list)  Internal Lab Sample Transfer Records and Tracking Sheets (describe or list)	S	_			<del></del> . <del></del>	
Airbills (No. of shipments) Chain-of-Custody Records Sample Tags Sample Log-In Sheet (Lab & DC1) Miscellaneous Shipping/Receiving Records (describe or list)  Internal Lab Sample Transfer Records and Tracking Sheets (describe or list)  Other Records (describe or list)		ll instrument ou	tput, including strip			
Airbills (No. of shipments) Chain-of-Custody Records Sample Tags Sample Log-In Sheet (Lab & DC1) Miscellaneous Shipping/Receiving Records (describe or list)  Internal Lab Sample Transfer Records and Tracking Sheets (describe or list)  Other Records (describe or list)					•,	
Airbills (No. of shipments) Chain-of-Custody Records Sample Tags Sample Log-In Sheet (Lab & DC1) Miscellaneous Shipping/Receiving Records (describe or list)  Internal Lab Sample Transfer Records and Tracking Sheets (describe or list)  Other Records (describe or list)						
Airbills (No. of shipments) Chain-of-Custody Records Sample Tags Sample Log-In Sheet (Lab & DC1) Miscellaneous Shipping/Receiving Records (describe or list)  Internal Lab Sample Transfer Records and Tracking Sheets (describe or list)  O. Other Records (describe or list)	. <u>ef</u>	PA Shipping/Rece	piving Documents		•	
Sample Log-In Sheet (Lab & DC1)  Miscellaneous Shipping/Receiving Records				<del></del> -		
Sample Log-In Sheet (Lab & DC1)  Miscellaneous Shipping/Receiving Records    (describe or list)  Internal Lab Sample Transfer Records and Tracking Sheets (describe or list)  Other Records (describe or list)		<del>-</del>	Records			
Internal Lab Sample Transfer Records and Tracking Sheets (describe or list)  Other Records (describe or list)	Sa	ample Log-In She				
Sheets (describe or list)  Other Records (describe or list)	Mi			S		
Sheets (describe or list)  Other Records (describe or list)			<del></del>			
Sheets (describe or list)  Other Records (describe or list)			· · · · · · · · · · · · · · · · · · ·	·		
				1 Tracking		
		<del></del>	<del></del>			
	 o. <u>c</u>	Other Records (d	escribe or list)			
	_		<del></del>			

CASE NO.	SDG	NO	SDG NOS. TO FOL	LLOW	
		SAS NO.			
11. Comments:					
			<del></del>		
•		•			,
Completed by:				·	
CLP Lab)	(Signature)		(Printed Name/	Title)	(Date)
erified by: _					
(CLP Lab)	(Signature)		(Printed Name/	Title)	(Date)
Audited by: _			<del></del>	· · · · · · · · · · · · · · · · · · ·	
(EPA)	(Signature)		(Printed Name/	Title)	(Date)

#### COVER PAGE - INORGANIC ANALYSES DATA PACKAGE

Lab	Name:			Contract:	<del></del>
Lab	Code: _		Case No.:	SAS No.:	_ SDG No.:
SOW	No.: _				
		EPA	Sample No.	Lab Sampl	e ID.
					<del></del>
					<del></del>
			<del></del>		<del></del>
					<del></del>
					<del></del>
		<del></del>			<del>-</del>
					<del></del>
			<del>,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,</del>		<del></del>
			<del></del>		<del></del>
			<del></del>	<del></del>	
		<del></del>			<del></del>
				<del></del>	<del></del>
		<del></del>		<del></del>	<del></del>
Were	ICP in	terele	ment corrections ap	oplied?	Yes/No
Were			nd corrections appl		Yes/No
	If yes	-were	raw data generated	before	-
	applic	ation o	of background corre	ections?	Yes/No
Comme	ents:				
					<del></del>
				<u>,, ,, ,, ,, ,, ,, ,, ,, ,, ,, ,, ,, ,, </u>	<del></del>
I cer	tify t	hat thi	is data package is	in compliance with	the terms and
condi than	the co	or tne nditior	contract, both ted as detailed above.	nnically and for c Release of the da	ompleteness, for oth ta contained in this
hardo	copy da	ta pac	tage and in the con	puter-readable dat	a submitted on diske
has b	een au	thorize	ed by the Laborator ollowing signature.	y Manager or the M	anager's designee, a
Signa	ture:_			Name:	
_					

U.S. EPA - CLP EPA SAMPLE NO. INORGANIC ANALYSIS DATA SHEET Lab Name: Contract: Lab Code: \_\_\_\_ SAS No.: \_\_\_ SDG No.: \_\_\_ Matrix (soil/water): \_\_\_\_ Lab Sample ID: Date Received: Level (low/med): % Solids: Concentration Units (ug/L or mg/kg dry weight): CAS No. | Analyte Concentration C M Q 7429-90-5 Aluminum 7440-36-0 Antimony\_ 7440-38-2 Arsenic\_\_ 7440-39-3 Barium 7440-41-7 Beryllium 7440-43-9 Cadmium 7440-70-2 Calcium 7440-47-3 Chromium\_ 7440-48-4 Cobalt\_\_\_\_ 7440-50-8 Copper 7439-89-6|Iron 7439-92-1 Lead 7439-95-4 Magnesium 7439-96-5 Manganese 7439-97-6 Mercury 7440-02-0 Nickel 7440-09-7 Potassium 7782-49-2 | Selenium\_ 7440-22-4 | Silver\_\_\_ 7440-23-5 Sodium 7440-28-0 Thallium 7440-62-2 | Vanadium |

	l					
Color	Before:		Clarity Before	<b>:</b>	Texture:	
Color	After:		Clarity After	***************************************	Artifacts:	
Comme	nts:		·		•	
-						
•		· · · · · · · · · · · · · · · · · · ·	·····	<del></del>		

7440-66-6 Zinc Cyanide

# 2A INITIAL AND CONTINUING CALIBRATION VERIFICATION

			<b></b>	<del></del>	_			
Lab Code: _	Ca	se No.:		SA	s No.:	<del></del>	\$	SDG No.:
Initial Cal	libration	Source:						
Continuing	Calibrati	on Sour	ce:	<del></del>				
			Concent	ration 1	Units: ug	/L		
	Initial	Calibr	ation	,	Continuin	g Calik	ration	
Analyte	True	Found	<b>%</b> R(1)	True	Found			%R(1)
Aluminum	<sub>1</sub> -	<del></del>	ı			11		-,
Antimony_		<del></del>						-
Arsenic		<del></del>						-
Barium —						[		-
Beryllium				•				-
Cadmium						1		-
Calcium		<del></del>						-
Chronium_		<del></del>				1		-
Cobalt		<del></del>				11		-
Copper				<del></del>				-
Iron		<del></del>						-
Lead								-
Magnesium								-
Manganese	<del></del>							-
Mercury	<del></del>			***************************************				-
Nickel								-
Potassium	<del></del>							-
Selenium								-1
Silver		<del></del>						-
	<del></del>							-
Sodium	ŀ							

(1) Control Limits: Mercury 80-120; Other Metals 90-110; Cyanide 85-115

#### 2B CRDL STANDARD FOR AA AND ICP

Lab Name:		Contract:	<del></del>	
Lab Code:	Case No.:	SAS No.:	SDG No.:	_
AA CRDL Standard S	ource:			
ICP CPDL Standard	Source:			

Concentration Units: ug/L

CKDL S	tandard fo	or AA		RDL Standa itial		Final	
True	Found	<b>≹</b> R	True	Found	<b>%</b> R	Found	%F
<u> </u>		,		ı ———	1	ı ————	,
·				l		l ———	
		<del></del>					
·		[———					
	<del></del>					<del></del>	
·				·		<del></del>	
							-
			<u> </u>	İ			
							<b> </b>
			1				
			<u> </u>				<del></del>
[.				·			
	}			· ———			
	)						
				True Found &R True		True Found &R True Found &R	True Found &R True Found &R Found

Control Limits: no limits have been established by EPA at this time

#### 3 BLANKS

	Lab Name:		Contract:	<del></del>
_	Lab Code:	Case No.:	SAS No.:	SDG No.:
	Preparation Blank	Matrix (soil/water): _	· ·	
	Preparation Blank	Concentration Units ()	ng/L or mg/kg):	

Analyte	Initial Calib. Blank (ug/L)	С		ti B C	nuing Cali lank (ug/L 2	br ) C		С	Prepa- ration Blank	С	M
Aluminum		.,-	]	, –	1	,-	1	-,-		,-	
Antimony_		· -		-		-	<del></del>	-1-1		-	
Arsenic	<del></del>	-		-	ļ ———	-1-		- -		1-1	I
Barium	<del></del>	-1-		1-	ļ <del></del>	-		-1-1		-	]—
Beryllium		1-	<del></del>	-	[ <del></del>	-		-1-1		1-1	1—
Cadmium		1-		1-		-		- -		-	
Calcium		-		-	ļ <del></del>	-	<del></del>	- -		-	]-
Chromium		-		-	]	-		-1-1		1-1	
Cobalt	<del></del>	<b> -</b>		-	ļ <del></del>	-	ļ ————	-1-1	ļ	-	-
Copper		[-]		-	l <del></del>	1-	l	-1-1		1-1	<b> -</b>
Iron		[ <b>-</b> ]		-	ļ <del></del>	-	<u> </u>	-!-!		-	
Lead		[-]		-		-	ļ <del></del>			1-1	
Magnesium		-		-		-	<b> </b>	- -		-	
Magnesium		[-]		-		-		- -		1-1	
		-		]_		1-	<del></del>	-1-1	\ <u></u>	-	
Mercury	<del></del>	1-1		-		-	<del></del>	-1-1	ļ <del></del>	-	
otassium		-		-		-		-1-1		-	
	<del></del>	-		-		-	<del></del>	-1-1		-	
Selenium_		-		<b>]</b> —		-		-1-1		1-1	<b> </b>
		1-1		<b> </b>		<b> </b>		-1-1			<b> </b>
Sodium		-		_	·	-		- -		-	
Thallium		-		-		-	·	- -			
/anadium_				<b> </b>		-		-		-	
Zinc		-		_		-		- -		_	_
Cyanide		_		_		_		- -		_	
		I I		l		l		_1_1	1	1	1

### ICP INTERFERENCE CHECK SAMPLE

Lab	Name:		Contract		-	
Lab	Code:	Case No.:	SAS No.:		SDG No.:	
ICP	ID Number:		ICS Source	ce:		

Concentration Units: ug/L

	T	rue		itial Fou	nd		al Found	·
	Sol.	Sol.	Sol.	Sol.		Sol.	Sol.	
Analyte	A	AB	A	AB	%R	A	AB	&R
Aluminum			<u> </u>		1			1
Antinony_								
Arsenic							}	
Barium								
Beryllium								
Cadmium		l l		<u> </u>	<b> </b>			
Calcium_		]				·	l	l
Chromium_						<u> </u>	·	<b> </b>
Cobalt					·		\	l
Copper		]	<b></b>					<b> </b>
Iron				]	<b> </b>	<b> </b>		<b> </b>
Magnesium							[ <del></del>	
Manganese							l	
Mercury								
Nickel -					ļ			
Potassium								
Selenium							·	
Silver				ļ ———				
Sodium			·				· · · · · · · · · · · · · · · · · · ·	
Thallium								
Vanadium								
Zinc								

# 5 A

		SPI	KE :	5A SAMPLE RECOVE	RY			EPA S		
Lab Name:	····		. <del></del>	Contract	::			]		
Lab Code:	b Code: Case		: _	SAS No.:				SDG No	.:	_
Matrix (so:	il/water)	:				1	<b>Le</b> vel	(low/me	: (£	_
k Solids fo	or Sample	:								
ı ————————————————————————————————————	Concer	tration Units	(uḍ	g/L or mg/kg	dry	weight	): <u> </u>	<del></del>	_	,
Analyte	Control Limit %R	Spiked Sampl Result (SSR	e ) C	Sample Result (SR)	c	Spi Added	ke (SA)	&R		M
_			-, -		_,_	\ <u></u>				_
Aluminum			_   _		- -	ļ			.   _	ļ
Antinony_			_ _		-1-				. _	<b> </b>
Arsenic			- -		-1-	<b></b>			.]_,	
Barium Beryllium			- -		- -				-   -	<b> </b>
Cadmium		<del></del>	- -	<del></del>	- -				-   -	-
Calcium			-1-1		-]-	J			- [ [	-
Chromium_			- -		- -	]			-	-
Cobalt			-]-		- -	]			-	-
Copper		<del></del>	- -	·····	-)-	J			-	-
Iron		<del></del>	- -	<del></del>	-]-	<u> </u>	<del></del>		-	
Lead	<del></del> }		_ _		- -	]			-	-
Magnesium	<del></del> ].		, ,	·	- -				-	-
Manganese	<del></del> ]·		- -		- -	l			1-1	-
Mercury	<del></del> j·		- -		- -				-	-
Nickel		······································	-1-1		- -	]			-	
Potassium	<del></del> [·		- -		- -		<del></del> [		1-1	_
Selenium_	<del></del> [.	<del></del>	-1-		- -		[		1-1	
Silver	[·		-1-1	<del></del>	- -				-	_
Sodium		<del></del>	- -		-1-			<del></del>	1-1	
Thallium			- -	· · · · · · · · · · · · · · · · · · ·	_ _		<del></del>		1-1	
Vanadium_			- -		_ _					_
Zinc										_
Cyanide									1_1	
			_ _		_[_					
•					_					
omments:										

# 5B

ab Name:			Contract:					-
•			•			·		
ab Code:	<del></del>	Case No.:	SAS No	• :		SDG No	•	
atrix (so	il/water)	):			Level	(low/me	d):	
					•			
		•						
		Concentra	tion Units: ug	/L		•		
	Control	V					T	
	Limit	Spiked Sample	Sample		Spike	İ	ł	l
Analyte	₹R	Result (SSR) (	Result (SR)	С		₹R	Q	М
	lI			. —			. _	
Aluminum_ Antimony_				-		<del></del>	- -	
Arsenic	<del></del>		.	-		<del></del>	-	
Barium				-			-	
eryllium			·   <del></del>	-			-	
admium	<del></del>		·	-		<del></del>	-	<b> </b> —
Calcium				-			-	
hromium				-			·[-	
Cobalt			1	-			-	
copper	<del></del>		1	-			1-	
ron			1	-			-	-
ead			i	-			-	_
Magnesium				-			-	
langanese								
ercury			ł.	-			-	
lickel —		1	1					
otassium				_			-	
elenium_								
Silver								
odium			1					
hallium_				I_I			-	
anadium				-				
inc				-			-	_
yanide				_				
- Name — I			1	ı — I			I = I	· —

6	
DUPLICATES	

-		וטס	PLICATES	EPA SAMPLE NO.
-	Lab Name:		Contract:	
	Lab Code:	Case No.:	SAS No.:	SDG No.:
	Matrix (soil/water):		Lev	rel (low/med):
-	% Solids for Sample:		% Solids for	Duplicate:

Concentration Units (ug/L or mg/kg dry weight): \_\_\_\_

Analyte	Control Limit	Sample (S)	С	Duplicate (D)	С	RPD	Q	M
Aluminum				\ <del></del>	- <sub>1</sub> -		-	1-
Antimony_		<del></del>	- -		- -		-	-
Arsenic -		<del></del>	- -		- -		-	-
Barium —		<del></del>	-1-1		-1-1		1-	1-
Beryllium		<del>· · · · · · · · · · · · · · · · · · · </del>	- -		- -		-	-
Cadmium			- -		- -		-	-
Calcium		<del></del>	- -		- -		-	-
Chromium			-1-1		_ _		-	-
Cobalt			_ -				_	
Copper			- -		- -		-	
Iron	<del></del>		_ _				1	
Lead					_/_/		1	
lagnesium			_ _		_(_(			1.
langanese			_ _				1	I
fercury					_ _			
Nickel					_ _			łΞ
Potassium			<u> </u>		_ _			ΙI
Selenium_					_ _			
Silver			_ _		_ _		1_	_
Sodium					_ _		_	_
Thallium_					_ _		1_	_
/anadium_			_ _		_ _		_	
Zinc			_ _		_ _			_
yanide					_   _			-

### 7 LABORATORY CONTROL SAMPLE

Lab Name:		Contract:		
Lab Code:	Case No.:	SAS No.:	SDG No.:	•
Solid LCS Source:	<del></del>			
Aqueous LCS Source:				

	Aqueous (ug/L)								
Analyte	True	Found	*R	True	Found	C	: Li	mits	*
Aluminum	·——	<del></del>		!		7 -	<u> </u>	1	-1
Antimony						-			
Arsenic			I.			-			-
Barium -			I			-			1
Beryllium_		(				-			'
Cadmium -						-			·
Calcium						-			-
Chronium						-			-
Cobalt						-			
Coppeir	<del></del>  -	·	-	<del></del>		-			-
Iron			——I·			-	ļ <del></del>		1-
Lead			·			-		<del></del>	·
Magnesium		<del></del>  -	-			-			·
fanganese_	<del></del> [-	<del></del> [-	(·			-		[ <del> </del>	·
fercury						-			·
lickel		<del></del>	-		<del></del>	-			·
otassium			-		·	-			
elenium		·	-			-			·
ilver						-			·
odium	-		———]·			-			· -
	<del></del>  -		—— -			-		<del></del>	.
hallium_	-					-			·
anadium			[-			_			<b> -</b>
inc	_		-			_			<b> </b> —
yanide	_					_			<b> </b>

# STANDARD ADDITION RESULTS

-	Lab Name:			Contract:		
	Lab Code:	<del></del>	Case No.:	SAS No.:	SDG	No.:

Concentration Units: ug/L

EPA	T										$\overline{}$
Sample No.	An	0 ADD ABS	1 A CON	DD ABS	2 A CON	DD ABS	3 A CON	DD ABS	Final Conc.	r	Q
											: -
											:  <u>-</u>
	-									<u></u>	- -
											- -
											-
											-
											:   _
											-
											-
	_										-
											-
											<u> -</u>
											-
											-
											-
		]									
	I.				<u></u>			l			1_

9
ICP SERIAL DILUTIONS

Lab Name:		Contract:	
Lab Code:	Case No.:	SAS No.:	SDG No.:
Matrim (soil/wa	ter):	Level	(low/med):

Concentration Units: ug/L

			Serial		<b>3</b>	T	
1	Initial Sample		Dilution	1	Differ-		]
Analyte	Result (I)	c	Result (S)	c	ence	Q	M
- 1	1	Ì	1	i			l
Aluminum		-		-1-1		-	1-
Antimony_		-		- -		-	
Arsenic -		-		- -		-	
Barium —		-		- -		-	_
Beryllium		-		-1-1		1-	-
Cadmium		-		- -		-	
Calcium		-	<del></del>	- -		-	_
Chromium		-		-		-	
Cobalt		-1	<del></del>	- -		-	_
Copper		-		- -		-	-
Iron		-	<del></del>	- -		-	
Lead		-1	<del></del>	- -		-	
Magnesium		-1	<del></del>	-   -			
Manganese		-[	[ <del></del>	-1-1		-	
Mercury		-)	<del></del>	- -		-	
Nickel -		-		- -		-	-
Potassium		-		- -		-	-
Selenium	<del></del>	-	[ <del></del>	-		-	-
Silver		-		1-1		1-1	-
Sodium		-[		- -		-	
Thallium		- [		-1-1	II	1-1	_
Vanadium_		-		- -		1-1	_
Zinc Zinc	[ <del></del> ]	-		- -		1-1	
Zine		-1		. -		1-1	

EPA SAMPLE NO.

# 10 INSTRUMENT DETECTION LIMITS (QUARTERLY)

<del></del>	Case No.:		<del></del>	3DG 1		
CP ID Number:		Dat	ie:	· · · · · · · · · · · · · · · · · · ·	-	
lame AA ID Nu	mber:					
urnace AA ID	Number:					
	)	Wave-	71	CDDI	TDT	
	3 - 3 - 4 -	length		CRDL	IDL	1,,
	Analyte	$(n\pi)$	ground	(ug/L)	(ug/L)	M
	Aluminum	<del></del>		200		-
	Antimony_			60		-
	Arsenic_			10		
	Barium			200		
	Beryllium			5_		. _
	Cadmium_			5		.
	Calcium			5000		-
	Chromium_		·	10		-
	Cobalt		-	50		·]—
	Copper		·	25 100	<del></del>	-
	Lead		·	3		-
	Magnesium			5000		-
	Manganese			15_		-
	Mercury			0.2		
	Nickel			40		
	Potassium			5000		.
	Selenium_			5_		.
	Silver			10	····	.]
	Sodium	<del></del>	-	5000		·
	Thallium			10 50		· —
	Zinc			20		-
	Cyanide			10		1
	<del>X   X   X   X   X   X   X   X   X   X</del>				<del></del>	1

# 11A ICP INTERELEMENT CORRECTION FACTORS (ANNUALLY)

	Ab Code: Case No.:			Contract: SDG No.:  Date:				
	*		· · · · · · · · · · · · · · · · · · ·		· · · · · · · · · · · · · · · · · · ·			
	Wave-		Interelement	Correction	Factors fo	or:		
Analyte	length (nm)	Al	Ca	Fe	Mg			
Aluminum Antimony Arsenic Barium Beryllium Cadmium Calcium Chromium Cobalt Copper Iron Lead Magnesium Manganese Mercury Nickel Potassium Selenium Silver Sodium Thallium Vanadium Zinc								

# 11B ICP INTERELEMENT CORRECTION FACTORS (ANNUALLY)

Lab Chue	<del></del>	Case No.: _	<del></del>	SAS No.:		·
ICP ID Numb	er:			Date:		
	Wave- length		Interelement	Correction	Factors fo	r:
Analyte	(nm)			<del></del>		
Aluminum_		<del> </del>	-1	ı	i <del></del>	
ADTIMODY					<del></del>	
Arsenic						
Barium			•			
Beryllium						
Cadmium						
Calcium						
Chromium						
Copait						
Copper					l	]
Iron						
Lead		<u> </u>	-	l		
Magnesium						
Manganese		<del></del>	.			
Mercury	]		.	J	<del></del>	]
Nickel						
Potassium						
Selenium_						
Silver		<del></del>				
mballine.	<del></del>		·			
Thallium Vanadium Zinc	<del></del> [					
7ino			·			
Zinc						<del></del> -
		l	.			l
omments:						

# 12 ICP LINEAR RANGES (QUARTERLY)

P ID Number:		1	Date:	
	Analyte  Aluminum_ Antimony_ Arsenic_ Barium Beryllium Cadmium_ Calcium_ Chromium_ Cobalt_ Copper_ Iron_ Lead_ Magnesium Manganese Mercury_ Nickel_ Potassium Selenium_ Silver_ Sodium_ Thallium_ Vanadium	Integ. Time (Sec.)	Concentration (ug/L)	M
	Zinc			

#### 13 PREPARATION LOG

-	Lab Name:	·	Contract:	_
	Lab Code:	Case No.:	SAS No.:	SDG No.:
•	Method:			

EPA	T	<del></del>	T
Sample No.	Preparation Date	Weight (gram)	Volume
No.	Date	(gram)	(mL)
		(32,	,
	·   <del></del>		
	·		
	·   <del></del>		
		·	
J			
	<del></del>		
			i
l			
	\ <u></u>		J
		1	
	<b>—————</b>		
		<del></del>	
	J	<del></del>	
İ		<u> </u>	
Ì			l
<del></del>			
<del></del>	<del></del>	<del></del>	<del></del>
<del></del>		<b></b>	
		[	
	<del></del>	·	

#### 14 ANALYSIS RUN LOG

Lab Name:											C	on	tr	ac'	t:	_				_							
Lab Code:			Case	Ио	.:	_			_		Ş	AS	N	٥.	: _					SI	DG	No	<b>5.</b> 1	: .			
Instrumen	t ID Numb	er: _	<del> –</del>								M	et!	ho	d:	_	_											
Start Dat	e:	<del></del>									E	nd	D	ate	e :	_				_							
EPA				Ī				-					A	na.	ly	tes	5					~					
	D/F	Time	* R	Ā	S	AS	В	BE	CD	C	CR	C	C	FE	P B	M G	M	H G	N I	K	S	A G	N A	T	V	Z N	
		.	ļ		_	_	_	_	_	_	_	_	_	_	_	_	_	_	_	_	_	_	_	_	_	_	_
			ļ <del></del>	-	-	-	-	-	-	_	<b>-</b>	_	_	-	-  -	-	-	-	<u>-</u>	_	_	-	-	-	_	-	-
				-	-	_	_	_	_	_	_	_	_	_	_	-	1 1	1 1	_	-	_	_	_	-	-	-	_
		·		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
					_	_	_	_	_	_	-	_	_	_	_	-	_	_	_	_	_		_	_	_		_
		·		-	-	-	-	-	-	-	-	-	-	-	-	-	1	-	_	1 1	-	-	-	-	-	-	-
				-		_	_	_	<u> </u>	_	1 1	_	_	_	1	_	1 1	- 1	_	- 1	-	_	_	_	_		1
		.		_	_	-	_	-	-	-	-	_	_	-	_		-	-	_	_	_	-	_	_	_	_	-
	I———	·		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
				_	_	_	_	_	+	_	1	_	_	1 1	_	_	_		1 1	_	-	_	_	_	_		
		.	ļ	-	<u> </u>	-	_	_	_	_	_		-	-		_	-	-	_		_	_	_	-	_	-	_
			<b></b>	[-	-	-	-	-	-	-	-	-	-	-	_	-	-	-	-	-	-	-	_	-	-	-	-
					_	_	1 1	-	_	_		-	_	_					_		_	_					
	<u> </u>		<b> </b> -	-	-	-	_	_	-	_	-	_	_	_	-	-	-	-	-	-	_	_	-	-	-	-	-
				-	-	-	-	-	-	-		_	_	_	_			-	-	-	-	-	-	-	-	-	-
					1	_	_		_			_	_	_			- 1								_		
<del></del>				-		-	-	_	_	_	_	-	_	_	_	_	_	_	_	_	-	_	_	-	-	-	-
	ļ <del></del>	ļ		-	-	-	-	-	-	-	-1	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
											_		_	_							_						
				-	_	_	_	_	_	_	_	_	_	_	_	_	ı	_	_	_	_	_	_	_	_	_	_
	ļ <del></del>			[-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-1
				-	-	-	-	-	-	-	-	-	-	-	-	- 1	- 1	-	-	-	-	-1	-	-	-	-	-
				_	_													_									
				_	_	_	_	_	_	_	_	_	_	_	_	_	_	_	_	_	_	_	-1	_	_	_	_
		<b> </b>		-	-	-	-	-1	-	-	-	-	-	-	-	-	-	-	-	-1	-	-	-	-	- -	- -	-
		1		ا ـــ ا	_	-1	_	_	_	1	1	_ 1	_	_ '	- 1	I	1	<b>—</b> I.	1	<b></b> - [	1	_1	!	!	<b></b> i.	_1.	_ 1

ATTACHMENT B

FIELD DATA SHEETS

PROJECT NO.		HULL					<u></u>	OF
		DAIL	Y FI	ELD	REPO	RT		
PROJECT		<del></del>			<del></del>	DATE_		
LOCATION					<del></del>	WEATH	ER	
CONTRACTOR	₹				<del></del>	TIME O	N-SITE	FROM
•								то
VISITORS ON	SITE				<del> </del>			· ·
				<u> </u>				
DESCRIPTION	of WORK							
					<del> ,</del>			
				<del></del>				
<del></del>					<del></del>	<del></del>		
								<del></del>
		<u> </u>						
				-				<del></del> -
		· · · · · ·			<u></u>			
								Hou
				1	UPMEN"	r USED.		Hou
				1	JIPMEN"	r USED.		Hou
				1	JIPMEN"	T USED.		Hou
				1	JIPMEN"	T USED.		Hou
				1	JIPMEN"	r USED.		HOU
				EQU	JIPMEN'	T USED.		Hou
				EQU		T USED.		HOU
				EQU		T USED.		HOU
				EQU		T USED.		HOU
				EQU		T USED.		HOU
				EQU		T USED.		HOU
				EQU		T USED.		HOU
				EQU		T USED.		HOU
				LAB	OOR		BY	

Hull & Associates, Inc.
6130 Wilcox Road
Dublin, Ohio 43016

#### FIELD DATA SHEET GROUND-WATER MONITORING WELL SAMPLING

			Well I.D	
Client:	·	Site Location:		
	Project No			
Air Temperature:	We	ather Conditions:	<u> </u>	
Type of Well Construction_				
Condition of Well circle (G	Good / Poor) if poor, specify			
Cap Locked (Yes / No)	Lock No.:			
	feet Total Depth			
Free Product circle (Yes / N	(o) Depth to Product			feet
roduct Thickness	feet			
ample Date	Sample No.:	<del></del>	1	<u>_</u>
PURGE		- I	G/WELL VOLUME	
VOLUME PURGED		EMP: pH	CONDUCTIVIT units	Y
(GALEONS)/WELL		AT:25°C		
VOLUME			AT 25°C	
NA	STATIC			
	1			
	2			
			ن	
	3	·		
				·
	4			
	5		•	
		G "		
Well Volume Equals:		Gallons	·	
	•	·		
ments:				
	Purge Water			

sbect	 oſ	
PDCCI	 o.	

#### HULL & ASSOCIATES, INC. WELL DEVELOPMENT FORM

Job Number Site: Developers:			Well No and Type:  Initial Total Depth (ft TOC):  Final Total Depth (ft TOC):  Weather:							at hrs. at hrs.
Date	Time	Purge Method	Pumping Rate	Volume Purged <sup>d</sup>	DTW•	pHt	Temp.	Spec. Cond. <sup>h</sup>	Turbidity <sup>t</sup> .	Comments
		~							-	
		<del></del>								
					<u> </u>			,		
	-							· · · · ·		
					-					
	-	-			-		-			
					1					
					-	-			,	

ľλ	111	$\cap$	FFF.	CO	ME	TI	<b>TZ</b>
111	***	v		$\sim$	1112		برورز

Start Time:

Volume Water:

End Time:

Volume Sediment:

- a. Top of casing.
- b. NAPL nonaqueous phase liquid.
- c. Gallons per minute.

- d. Cumulative gallons,
- c. Depth to water.

- f. Standard units.
- g. \*C, unless \*F noted.

- h. Specific conductance, µmhos/cm (or µS/cm).
- i. Visual unless otherwise noted,

# Water Level and Interface Measurement Sheet

HAI Project #	•	Site Location
HAI Site Personnel		Date

Well Name	Static Water Level	Depth to Bottom	Volume	Depth lo Product	Product Thickness	Comments
					·	
·						
	·					
		<u> </u>				
					,	
		•				

HULL & ASSOCIATES, INC. SOIL BORING FIELD LOG

												· · · · · · ·		
LOC/	ATION OF BO	RING:							PROJ No:	CLIENT:		LOCATION:	·	···
									PROJECT:					BORING No:
									DRILLING METHOD:				· · · · · · · · · · · · · · · · · · ·	]
									SAMPLING METHOD:				<del></del>	SHEET
										LCVCBCUNC				1
									PID/FID CALIBRATION,BA	CAGROUND.				OF OF
<u>.</u>									WEATHER:		r <del></del>	1	1	LING
٩T	RACTOR:								WATER LEVEL FROM:			<u> </u>	START	FINISH
G	ED BY:				DAT	E:			TIME:				TIME:	TIME:
ECI	KED BY:				DAT	E:			DATE;				DATE:	DATE:
TUI					ELE	VATIO	N;	•					1	
	DRV/REC.	SMPL No /DEPTH	NOW COUNTS	PIO/FID (ppm)	DEFTH W FEET	SAUPLE	7 -	TIBM	DEPTH: NOTES (SURFACE COND	itión, lab sóil s	NUPLE NUMBERS	SOIL DRUMS, I	TC.):	•
	ļ		-			-	-							
_		<del></del>	-	-		+	$\dashv$							<del></del>
			-{			-	7						<del></del>	
1			1			1	1							
J			]											
					$\top$									
ſ					<u>                                     </u>		_				•			
			}				1							
l			<u> </u>			1_	1							·
			1		1	_	.			·····				
_			<u> </u>		J	$\perp$	1		<del> </del>			<u> </u>		
L						<u>_</u>								
_					_		1							
						L			· · · · · · · · · · · · · · · · · · ·					
				·	<u> </u>				·					
					1	Ш								
					T .									
-	,		1		]	П								
	—— -					$\vdash$		卜		<del></del>	<del></del>	<del></del>		
-	—— <u>-</u> -					$\dashv$		\ <u></u>		-				
-				<del></del>	<del> </del>	{	-	-			<del></del>			
_				į		_	}							
_								$\perp$						
			ł	1	[	$\Box$		L						
				ſ	<u> </u>									
_			1	ĺ	-									
_			-					-	- <del></del>					
_			-		}	-		-						
_					-			$\vdash$					<del> </del>	
				- 1	L	_		_					· · · · · · · · · · · · · · · · · · ·	
	·· .		_	l		_							<u></u>	
_			-[-			_								
_			1		Γ									
_						$\neg$								
-	<del></del>				 	$\dashv$		-						
_					-			-					·	
			-	1	_	$\dashv$		<u> </u>						
			_					_		<del></del>	<del></del>			
					_ [_	_								
				}		_								
	i	1	f		_(_									

# Hull & Associates, Inc.

□ Joledo

() Mason

□ Dublin

### CHAIN OF CUSTODY RECORD

☐ Warrensville Heights

PAGE\_\_\_\_ OF\_\_\_\_

NO. 1122

Phone: (514)793-8777 Phone:	nros Street Dhio 43606 (419)241-7171 (19)241-3117	4700 Duke Moson, Oh Phone: (5 FAX: (51	n Drive, Su No 45040 13)459–96 3)459–986	tta 172 4 77 1 9 1	1949 Galaxy Parkway, Warrensville Heights, Ol Phone: (218)514-7104 FAX: (216)514-7104	hlo 44128								ANALY	refe					
REPORT TO:						PRESE	RVATIV	FS	-/		_		7	ANALI	7	7	<del></del>			
Client: Site: Project#: Samplers:	Phase:			SAMPLE MATE 1, WATER 2, SOR, 3, SLUDGE 4, OR, 6, TISSUE 8, SEDIME OTHER:	A - Cool only, <# ( B - HNO3 pH<2 C - H2SO4 pH<2 D - NoOH pH>12	METALS  F - FILTE N - NOT B - BOTH									7					
SAMPLE IDENTIFICATION		SAMPLE MATRIX	NO. OF CONT.	METALS	SAMPLING DATE/TIME	/											COMMENTS			
					/															
				•																
			1				1													
	<del></del>		<del> </del>	1			1													
		<del>                                     </del>					+-	<del>                                     </del>		1	<u> </u>	1	1 -	<del>                                     </del>		1				
	<del></del>	1		1			1	†-		1	1	+-	1	1	1	1	<del> </del>			
		+	<del>                                     </del>	+		_	$\top$	1-		1	-	<del> </del>	+-		1	+-	<del> </del>			
RELINQUISHED BY:					RECEIVED BY:			DATE:					Deliver To:							
RELINQUISHED BY:	DATE:	DATE:			RECEIVED BY:			DATE:				Method of Delivery:								
RELINQUISHED BY:	DATE:		F	RECEIVED	DATE:				Airbill Number:  NOTES:											
COOLER TEMPERATURE AS RECEIVED 'C :	TIME:	DISTRIBUTION:  C WHITE — LAB USE (AUST YELLOW — LAB USE (AUST YELLOW — LAB USE (AUST NEW HE)																		