



US EPA RECORDS CENTER REGION 5



Headquarters/
Laboratory
705 front street
toledo, ohio 43605
phone: (419) 693-5307
fax: (419) 691-0418

Environmental
Laboratory
1632 enterprise parkway
twinsburg, ohio 44087
phone: (216) 425-8200
fax: (216) 425-1349

Sales Office
129 dixie way south
south bend, indiana 46637
phone: (219) 272-2200

Sales Office
77 west elmwood drive
suite 123
dayton, ohio 45459
phone: (513) 438-1711

Medical Products Division
615 front street
toledo, ohio 43605
phone: (419) 693-5307
fax: (419) 691-1227

October 25, 1996

**Ecology and Environment, Inc.
6777 Engle Road
Cleveland, OH 44130**

ATTN.: Mr. Dave Hendren

RE: QA II Documentation for KJ5100/S05-9610-805

Dear Mr. Hendren:

Enclosed is the QA II documentation requested for Project Name S05-9610-805 Project Number KJ5100 (B.E.C. Sample number 96T14933 through 96T14944).

If you have any further questions regarding the Quality Control data please do not hesitate to contact me at 1-419-693-5307.

Sincerely,

A handwritten signature in black ink, appearing to read "John F. Blair".

**John F. Blair
Quality Assurance Manager**

Enclosures



QA II Documentation

Project : S05-9610-805

Project Number: KJ5100

General Chemistry & Metals Analysis

Prepared for:
Ecology & Environment, Inc.
6777 Engle Road
Cleveland, OH 44130

ATTN: Mr. Dave Hendren

Submitted by:
B.E.C. Laboratories, Inc.
705 Front Street
Toledo, Ohio 43605

Table of Contents

| | |
|--|-----|
| Narrative | I |
| Chain of Custody | II |
| Analytical Results | III |
| Metals Raw Data & QC Reports | IV |
| Metals ICP Calibration Curves | V |
| General Chem. Raw Data & QC Reports | VI |

SECTION I

Narrative

Narrative
Ecology & Environment Inc.

Project: KJ5100 505-9610-805

BEC Sample Numbers 96T14933 to 96T14944

I. Chain of Custody

- 1) Samples were hand delivered in a cardboard box at ambient temperature.
- 2) Chain of Custody accompanied the samples and was properly filled out.
- 3) Chain of Custody seals were not used for this project.
- 4) Copies of the Chain of Custody can be found in section II of this document.
- 5) Copies of the analytical reports can be found in section III of this document.

II. pH Analysis

- 1) Due to the nature of these samples it was necessary to add water to determine the pH in all cases.
- 2) All pH's were either very high or very low.
- 3) pH meter calibration for 10/15/96 is included. Due to the low and high reading the meter was also checked @ a pH of 1.00 and 13.00 and found to be in calibration.

III. Cyanide & Reactive Cyanide & Flash Point

- 1) This analysis was run on two of the samples.
- 2) Batch Q.C. for total Cyanide showed a method blank of less than 0.1 mg/kg, method spike 82.2% recovery, duplicate results of 2.79 and 2.82 mg/kg for an RPD of 1.07% and a matrix spike recovery of 100.1%.
- 3) Batch Q.C. for Reactive Cyanide showed a method blank of less than 0.1 mg/kg, method spike 22.3% recovery, duplicate results of less than 0.1 mg/kg and a matrix spike recovery of 0.6%. Although these recoveries are low this is not unusual for Reactive Cyanide analysis.

4) Flash point was run on 96T14942 with a result of less than 74 degrees F p-xylene was analyzed on the instrument with a result of 79 degrees F which is within the 80 degree F +/- 1.5 degree F limits.

IV. Chloride Analysis

1) This was run on samples 96T14934 and 96T14936.

2) Batch Q.C. gave the following results: Method Blank less than 1 mg/L, method spike 108/4%, duplicates 19, 19 mg/l for an RPD of 0.0% and a matrix spike recovery of 107.1%.

3) Large dilutions were also necessary for this analyte due to high level of Chloride content.

V. Sulfate Analysis

1) This was run on sample number 96T14933.

2) Batch Q.C. gave the following results: Method Blank less than 2 mg/l, method spike 96.7%, duplicates 1720, 1740 mg/l for an RPD of 1.16% and a matrix spike of 95.5%.

3) As for other analytes dilution was necessary due to the high level of sulfate in the sample.

VI. Metals Analysis

1) Sample 96T14939 was run for Total Cr, Ni, and Zn. Dilution was necessary for Zn analysis.

2) ICP curves are included.

3) Batch Q.C. was as follows:

| Element | Method Blank | Method Spike | R1 mg/kg | R2 mg/kg | RPD | Matrix Spike |
|---------|--------------|--------------|----------|----------|-------|--------------|
| Cr | < 0.05 | 81.1% | 3.5 | 3.1 | 12.5% | 77.6% |
| Ni | < 0.05 | 92.3% | 5.8 | 4.8 | 18.9% | 71.6% |
| Zn | < 0.05 | 85.5% | 8.6 | 8.3 | 3.6% | 75.1% |

4) Printouts for the actual sample analysis are included.

Approved: John F. Blair Date: 10/25/96

Title: Qual. Assur. Manager



Section II

Chain of Custody

ENVIRONMENTAL PROTECTION AGENCY
Office of Enforcement

77 West Jackson Boulevard
Chicago, Illinois 60604

CHAIN OF CUSTODY RECORD

| PROJ. NO. | PROJECT NAME | CHAIN OF CUSTODY RECORD | | | | | | Activity Code: |
|---|--------------|-------------------------|------------|------|------------------|--|-------------|----------------|
| STA. NO. | DATE | TIME | COM. P. | GRAB | STATION LOCATION | NO. OF CON- TAINERS | TAG NUMBERS | MATRIX |
| D01 | 10/1/98 | 1000 | X | X | DRUM 01 - Room 1 | 1 | 937 | Liquid |
| D02 | 10/1/98 | 1010 | X | X | DRUM 02 - Room 1 | 1 | 934 | Liquid |
| D03 | 10/1/98 | 1012 | X | X | DRUM 03 - Room 1 | 1 | 935 | Liquid |
| D04 | 10/1/98 | 1015 | X | X | DRUM 04 - Room 1 | 1 | 936 | Liquid |
| V05 | 10/1/98 | 1018 | X | X | Vat 05 - Room 1 | 1 | 937 | Turnaround: |
| V06 | 10/1/98 | 1021 | X | X | Vat 06 - Room 1 | 1 | 938 | 14 day Verbal |
| V07 | 10/1/98 | 1025 | X | X | Vat 07 - Room 1 | 1 | 939 | 21 day hard |
| V08 | 10/1/98 | 1115 | X | X | Vat 08 - Room 2 | 1 | 940 | Copy |
| V09 | 10/1/98 | 1120 | X | X | Vat 09 - Room 2 | 1 | 941 | Liquid |
| D10 | 10/1/98 | 1130 | X | X | DRUM 10 - Room 3 | 1 | 942 | Liquid |
| D11 | 10/1/98 | 1145 | X | X | DRUM 11 - Room 3 | 1 | 943 | Solid |
| V12 | 10/1/98 | 1200 | X | X | VAT 12 - Room 3 | 1 | 944 | Liquid |
| Attn: Dave Hendren | | | | | | Ecology Environment, Inc. | | |
| Send Invoice to Dave Hendren | | | | | | 33 N. Dearborn Ste 900 | | |
| Consultant Cleve OH. Per David H. Hart | | | | | | 312-463-0791 Chicago IL 60602 | | |
| Relinquished by: (Signature) | | | | | | Ship To: | | |
|  | | | | | | BEC Lab 705 Front Street Toledo OH 43605 | | |
| Relinquished by: (Signature) | | | | | | ATTN: N. Huff (419-693-5307) | | |
| Relinquished by: (Signature) | | | | | | Airbill Number | | |
|  | | | | | | Date / Time | | |
| Received by: (Signature) | | | | | | Date / Time | | |
| Date / Time | | | | | | Received for Laboratory by: | | |
| (Signature) | | | | | | (Signature) | | |
| Date / Time | | | | | | Date / Time | | |
| 10/1/98 8:47 | | | | | | Received by: (Signature) | | |

Distribution: White - Accompanies Shipment; Pink - Coordinator Field Files; Yellow - Laboratory File

Chain of Custody Seal Numbers


5 11 E 1 C

SECTION III

Analytical Results

Ecology and Environment, Inc.
6777 Engle Road

Cleveland, OH 44130
ATTN: Dave Hendren



HEADQUARTERS/LABORATORY
705 FRONT STREET
TOLEDO, OHIO 43605
PHONE: (419) 693-5307
FAX: (419) 691-0418

ENVIRONMENTAL LABORATORY
1632 ENTERPRISE PARKWAY
TWINSBURG, OHIO 44087
PHONE: (216) 425-8200
FAX: (216) 425-1349

| | |
|----------|----------|
| lab no. | 96T14933 |
| p.o. no. | |

rev: 0

SAMPLE

DESCRIPTION: Project # KJ5100 - S05-9610-805 -- grab - liquid - Drum 01 - Room 1 -
Station # D01 - 10/11/96 @ 1000

RESULTS:

| <u>ANALYTE</u> | <u>METHOD</u> | <u>RESULT</u> |
|----------------|----------------|---------------|
| Sulfate | SW-846, 9038 | 820,000 mg/Kg |
| pH in 10% S.U. | EPA-600, 150.1 | less than 1 * |

* Due to the characteristics of the sample it was necessary to add water for pH determination.

date completed:
10/16/96 tech:
DWL

approved by:

Ecology and Environment, Inc.
6777 Engle Road

Cleveland, OH 44130
ATTN: Dave Hendren



HEADQUARTERS/LABORATORY
705 FRONT STREET
TOLEDO, OHIO 43605
PHONE: (419) 693-5307
FAX: (419) 691-0418

ENVIRONMENTAL LABORATORY
1632 ENTERPRISE PARKWAY
TWINSBURG, OHIO 44087
PHONE: (216) 425-8200
FAX: (216) 425-1349

| | |
|----------|----------|
| lab no. | 96T14934 |
| p.o. no. | |

rev: 0

SAMPLE

DESCRIPTION: Project # KJ5100 - S05-9610-805 - grab - liquid - Drum 02 - Room 1 -
Station # D02 - 10/11/96 @ 1010

RESULTS:

| <u>ANALYTE</u> | <u>METHOD</u> | <u>RESULT</u> |
|----------------|---------------|---------------|
| pH in 10% S.U. | SW-846, 9040 | less than 1 * |
| Total Chloride | SW-846, 9252 | 380,000 mg/L |


* Due to sample characteristics it was necessary to add water for
the pH determination.

date completed:
10/13/96

tech:
CPY/TW

approved by:

Ecology and Environment, Inc.
6777 Engle Road



Cleveland, OH 44130
ATTN: Dave Hendren

HEADQUARTERS/LABORATORY
705 FRONT STREET
TOLEDO, OHIO 43605
PHONE: (419) 693-5307
FAX: (419) 691-0418

ENVIRONMENTAL LABORATORY
1632 ENTERPRISE PARKWAY
TWINSBURG, OHIO 44087
PHONE: (216) 425-8200
FAX: (216) 425-1349

lab no.
96T14935
p.o. no.

rev: 0

SAMPLE

DESCRIPTION: Project # KJ5100 - S05-9610-805 - grab - liquid - Drum 03 - Room 1 -
Station # D03 - 10/11/96 @ 1012

RESULTS:

| <u>ANALYTE</u> | <u>METHOD</u> | <u>RESULT</u> |
|----------------|---------------|-------------------|
| pH in 10% S.U. | SW-846, 9040A | greater than 13 * |


* Due to the characteristics of the sample it was necessary to add water for
the pH determination.

| | | |
|-----------------------------|--------------|------------------|
| date completed: 10/15/96 | tech: CSV | approved by: |
|-----------------------------|--------------|------------------|

All reports are submitted as confidential communications. Authorization for duplication in whole or in part is reserved pending our written approval, as a mutual protection.

Ecology and Environment, Inc.
6777 Engle Road

Cleveland, OH 44130
ATTN: Dave Hendren



HEADQUARTERS/LABORATORY
705 FRONT STREET
TOLEDO, OHIO 43605
PHONE: (419) 693-5307
FAX: (419) 691-0418

ENVIRONMENTAL LABORATORY
1632 ENTERPRISE PARKWAY
TWINSBURG, OHIO 44087
PHONE: (216) 425-8200
FAX: (216) 425-1349

| | |
|----------|----------|
| lab no. | 96T14936 |
| p.o. no. | |

rev: 0


SAMPLE

DESCRIPTION: Project # KJ5100 - S05-9610-805 - grab - liquid - Drum 04 - Room 1 -
Station # D04 - 10/11/96 @ 1015


RESULTS:

| <u>ANALYTE</u> | <u>METHOD</u> | <u>RESULT</u> |
|----------------|----------------|---------------|
| pH in 10% S.U. | EPA-600, 150.1 | 11.67 * |
| Total Chloride | SW-846, 9252 | 170,000 mg/L |

* Due to sample characteristics it was necessary to add water for the pH determination.

| | | |
|-----------------------------|------------------|---|
| date completed: 10/15/96 | tech: CGY/CJE | approved by:  |
|-----------------------------|------------------|---|

Ecology and Environment, Inc.
6777 Engle Road



Cleveland, OH 44130
ATTN: Dave Hendren

HEADQUARTERS/LABORATORY
705 FRONT STREET
TOLEDO, OHIO 43605
PHONE: (419) 693-5307
FAX: (419) 691-0418

ENVIRONMENTAL LABORATORY
1632 ENTERPRISE PARKWAY
TWINSBURG, OHIO 44087
PHONE: (216) 425-8200
FAX: (216) 425-1349

| | |
|----------|----------|
| lab no. | 96T14937 |
| p.o. no. | |

rev: 0

SAMPLE

DESCRIPTION: Project # KJ5100 - S05-9610-805 - grab - liquid - Vat 05 - Room 1 -
Station # V05 - 10/11/96 @ 1018

RESULTS:

| <u>ANALYTE</u> | <u>METHOD</u> | <u>RESULT</u> |
|----------------|---------------|---------------|
| pH in 10% S.U. | SW-846, 9040A | 12.67 * |

* Due to the characteristics of the sample it was necessary to add water
for the pH determination.

| | | |
|-----------------------------|--------------|------------------|
| date completed: 10/15/96 | tech: DHH | approved by: |
|-----------------------------|--------------|------------------|

Ecology and Environment, Inc.
6777 Engle Road

Cleveland, OH 44130
ATTN: Dave Hendren



HEADQUARTERS/LABORATORY
705 FRONT STREET
TOLEDO, OHIO 43605
PHONE: (419) 693-5307
FAX: (419) 691-0418

ENVIRONMENTAL LABORATORY
1632 ENTERPRISE PARKWAY
TWINSBURG, OHIO 44087
PHONE: (216) 425-8200
FAX: (216) 425-1349

| | |
|----------|----------|
| lab no. | 96T14938 |
| p.o. no. | |

rev: 0

SAMPLE

DESCRIPTION: Project # KJ5100 - S05-9610-805 - grab - liquid - Vat 06 - Room 1 -
Station # V06 - 10/11/96 @ 1021

RESULTS:

| <u>ANALYTE</u> | <u>METHOD</u> | <u>RESULT</u> |
|----------------|---------------|---------------|
| pH in 10% S.U. | SW-846, 9040A | less than 1 * |

* Due to the characteristics of the sample it was necessary to add water for
the pH determination.

| | | |
|-----------------------------|-------------|---|
| date completed: 10/15/96 | tech: CV | approved by: A handwritten signature in black ink, appearing to read "Jim Huff". |
|-----------------------------|-------------|---|

Ecology and Environment, Inc.
6777 Engle Road



Cleveland, OH 44130
ATTN: Dave Hendren

HEADQUARTERS/LABORATORY
705 FRONT STREET
TOLEDO, OHIO 43605
PHONE: (419) 693-5307
FAX: (419) 691-0418

ENVIRONMENTAL LABORATORY
1632 ENTERPRISE PARKWAY
TWINSBURG, OHIO 44087
PHONE: (216) 425-8200
FAX: (216) 425-1349

| | |
|----------|----------|
| lab no. | 96T14939 |
| p.o. no. | |

rev: 0

SAMPLE

DESCRIPTION: Project # KJ5100 - S05-9610-805 - grab - solid - Vat 07 - Room 1 -
Station # V07 - 10/11/96 @ 1025

RESULTS:

| <u>ANALYTE</u> | <u>METHOD</u> | <u>RESULT</u> |
|----------------|---------------|---------------|
| Total Chromium | SW-846, 6010A | 49 mg/Kg |
| Total Nickel | SW-846, 6010A | 43 mg/Kg |
| Total Zinc | SW-846, 6010A | 10,000 mg/Kg |

| | | |
|-----------------|---------|--------------|
| date completed: | tech: | approved by: |
| 10/17/96 | KCB/ERS | |

All reports are submitted as confidential communications. Authorization for duplication in whole or in part is reserved pending our written approval, as a mutual protection

Ecology and Environment, Inc.
6777 Engle Road

Cleveland, OH 44130
ATTN: Dave Hendren



HEADQUARTERS/LABORATORY
705 FRONT STREET
TOLEDO, OHIO 43605
PHONE: (419) 693-5307
FAX: (419) 691-0418

ENVIRONMENTAL LABORATORY
1632 ENTERPRISE PARKWAY
TWINSBURG, OHIO 44087
PHONE: (216) 425-8200
FAX: (216) 425-1349

| | |
|----------|----------|
| lab no. | 96T14940 |
| p.o. no. | |


rev: 0

SAMPLE

DESCRIPTION: Project # KJ5100 - S05-9610-805 - grab - liquid - Vat 08 - Room 2 -
Station # V08 - 10/11/96 @ 1115

RESULTS:

| <u>ANALYTE</u> | <u>METHOD</u> | <u>RESULT</u> |
|------------------|-----------------|---------------------|
| Total Cyanide | SW-846, 9010 | 6500 mg/Kg |
| Reactive Cyanide | SW-846, 7.3.3.2 | less than 0.1 mg/Kg |

| | | |
|-----------------------------|------------------|---|
| date completed: 10/11/96 | tech: CSY/LMV | approved by:  |
|-----------------------------|------------------|---|

Ecology and Environment, Inc.
6777 Engle Road

Cleveland, OH 44130
ATTN: Dave Hendren



HEADQUARTERS/LABORATORY
705 FRONT STREET
TOLEDO, OHIO 43605
PHONE: (419) 693-5307
FAX: (419) 691-0418

ENVIRONMENTAL LABORATORY
1632 ENTERPRISE PARKWAY
TWINSBURG, OHIO 44087
PHONE: (216) 425-8200
FAX: (216) 425-1349

| | |
|----------|----------|
| lab no. | 96T14941 |
| p.o. no. | |

rev: 0

SAMPLE

DESCRIPTION: Project # KJ5100 - S05-9610-805 - grab - liquid - Vat 09 - Room 2 -
Station # V09 - 10/11/96 @ 1120

RESULTS:

| <u>ANALYTE</u> | <u>METHOD</u> | <u>RESULT</u> |
|----------------|---------------|---------------|
| pH in 10% S.U. | SW-846, 9040A | less than 1 * |

* Due to the characteristics of the sample it was necessary to add water for
the pH determination.

| | | |
|-----------------|-------|--------------|
| date completed: | tech: | approved by: |
| 10/13/96 | CSY | |

Ecology and Environment, Inc.
6777 Engle Road



Cleveland, OH 44130
ATTN: Dave Hendren

HEADQUARTERS/LABORATORY
705 FRONT STREET
TOLEDO, OHIO 43605
PHONE: (419) 693-5307
FAX: (419) 691-0418

ENVIRONMENTAL LABORATORY
1632 ENTERPRISE PARKWAY
TWINSBURG, OHIO 44087
PHONE: (216) 425-8200
FAX: (216) 425-1349

| | |
|----------|----------|
| lab no. | 96T14942 |
| p.o. no. | |

rev: 0

SAMPLE

DESCRIPTION: Project # KJ5100 - S05-9610-805 - grab - liquid - Drum 10 - Room 3 -
Station # D10 - 10/11/96 @ 1130

RESULTS:

| <u>ANALYTE</u> | <u>METHOD</u> | <u>RESULT</u> |
|-------------------------|---------------|----------------|
| Flash Point - SETAFLASH | SW-846, 1020 | less than 74°F |

date completed:
10/15/96

tech:
LMV

approved by:

A handwritten signature in black ink, appearing to read "Jim Neff".

Ecology and Environment, Inc.
6777 Engle Road



Cleveland, OH 44130
ATTN: Dave Hendren

HEADQUARTERS/LABORATORY
705 FRONT STREET
TOLEDO, OHIO 43605
PHONE: (419) 693-5307
FAX: (419) 691-0418

ENVIRONMENTAL LABORATORY
1632 ENTERPRISE PARKWAY
TWINSBURG, OHIO 44087
PHONE: (216) 425-8200
FAX: (216) 425-1349

| | |
|----------|----------|
| lab no. | 96T14943 |
| p.o. no. | |

rev: 0

SAMPLE

DESCRIPTION: Project # KJ5100 - S05-9610-805 - grab - solid - Drum 11 - Room 3 -
Station # D11 - 10/11/96 @ 1145

RESULTS:

| ANALYTE | METHOD | RESULT |
|------------------|-----------------|--------------|
| Total Cyanide | SW-846, 9010 | 53,000 mg/Kg |
| Reactive Cyanide | SW-846, 7.3.3.2 | 470 mg/Kg |

date completed:
10/16/96

tech:
CEN/LBN

approved by:

Ecology and Environment, Inc.
6777 Engle Road

Cleveland, OH 44130
ATTN: Dave Hendren



HEADQUARTERS/LABORATORY
705 FRONT STREET
TOLEDO, OHIO 43605
PHONE: (419) 693-5307
FAX: (419) 691-0418

ENVIRONMENTAL LABORATORY
1632 ENTERPRISE PARKWAY
TWINSBURG, OHIO 44087
PHONE: (216) 425-8200
FAX: (216) 425-1349

| | |
|----------|----------|
| lab no. | 96T14944 |
| p.o. no. | |

rev: 0


SAMPLE

DESCRIPTION: Project # KJ5100 - S05-9610-805 - grab - liquid - Vat 12 - Room 3 -
Station # V12 - 10/11/96 @ 1200

RESULTS:

| <u>ANALYTE</u> | <u>METHOD</u> | <u>RESULT</u> |
|----------------|---------------|---------------|
| pH in 10% S.U. | SW-846, 9040A | less than 1 * |

* Due to the characteristics of the sample it was necessary to add water for
the pH determination.

| | | |
|-----------------------------|--------------|---|
| date completed: 10/16/96 | tech: CSY | approved by:  |
|-----------------------------|--------------|---|

SECTION IV
Metals Raw Data & QC Data

Quality Control Data
Total Metals

| | | | |
|---------------------------|--|------------------------|------------------|
| Client: | <u>Ecology & Environment, Inc.</u> | Date Sampled : | <u>10/11/96*</u> |
| Project: | <u>S05-9610-805</u> | Date Received : | <u>10/15/96</u> |
| Project Number: | <u>KJ5100</u> | Date Analyzed : | <u>10/17/96</u> |
| Sample No.: | <u>Vat 07 – Room 1 – Station #V07</u> | | |
| B.E.C. Sample No.: | <u>96T14939</u> | | |
| Batch #: | <u>MET 76-41</u> | | |

| Metals | Units | Method | Method | Duplicate | | RPD | Matrix |
|---------------|--------------|---------------|---------------|------------------|-----|------------|---------------|
| | | Blank mg/l | Spike ** | R1 | R2 | | |
| Chromium | mg/kg | LT 0.05 | 81.1% | 3.5 | 3.1 | 12.5% | 77.6% |
| Nickel | mg/kg | LT 0.05 | 92.3% | 5.8 | 4.8 | 18.9% | 71.6% |
| Zinc | mg/kg | LT 0.05 | 85.5% | 8.6 | 8.3 | 3.6% | 75.1% |
| | | | | | | | |
| | | | | | | | |
| | | | | | | | |
| | | | | | | | |

* Date per Chain-of-Custody

LT – Less than

NA – RPD results are not applicable when duplicate results are below the detection limit.

** Method Spike is a lab characterized imitation soil.

Approved: J. F. Blair
Date: 10/25/96

E.C. Laboratories
7:17:43 16 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 31

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|--|-------|-------|--------|-------|-------|-------|-------|---|
| < Sample ID: 96T14939/3 10-16-96 Seq: 497 17:17:43 16 Oct 1996 ICP | | | | | | | | |
| g1 | .024 | PPM | .003 | .020 | .025 | .027 | .023 | |
| 2 | .031 | PPM | .047 | .055 | -.039 | .062 | .047 | |
| 2 | .001 | PPM | .000 | .001 | .001 | .000 | .001 | |
| a1 | .190 | PPM | .001 | .191 | .189 | .190 | .190 | |
| d1 | .041 | PPM | .002 | .044 | .040 | .040 | .040 | |
| 1 | .012 | PPM | .004 | .006 | .016 | .012 | .014 | |
| C 4 | 1.11 | PPM | .007 | 1.12 | 1.10 | 1.11 | 1.12 | |
| b1 | .229 | PPM | .012 | .219 | .233 | .245 | .221 | |
| 1 | .608 | PPM | .005 | .609 | .607 | .603 | .614 | |
| 1 | -.051 | PPM | .033 | -.027 | -.097 | -.054 | -.027 | |
| 11 | .029 | PPM | .039 | .079 | -.000 | .042 | -.003 | |
| N 3 | .975 | PPM | .011 | .979 | .974 | .961 | .987 | |
| 1 | .025 | PPM | .015 | .044 | .030 | .009 | .017 | |
| 2 | -.006 | PPM | .004 | -.011 | -.006 | -.003 | -.004 | |
| u1 | 32.3 | PPM | .237 | 32.3 | 32.1 | 32.2 | 32.7 | |
| 2 | 91.9 | PPM | .529 | 92.1 | 91.5 | 91.5 | 92.6 | |
| 1 | 191. | PPM | .907 | 192. | 190. | 191. | 192. | |
| 13 | 1.72 | PPM | .014 | 1.72 | 1.72 | 1.70 | 1.74 | |
| 2 | 3290 | PPM | 33.1 | 3340 | 3290 | 3260 | 3270 | |
| 3 | 23.9 | PPM | .151 | 23.9 | 23.8 | 23.7 | 24.1 | |
| g4 | 6.60 | PPM | .053 | 6.60 | 6.63 | 6.53 | 6.66 | |
| 1 | 6.38 | PPM | .448 | 6.28 | 5.99 | 6.21 | 7.02 | |
| 1 | .476 | PPM | .012 | .481 | .473 | .489 | .462 | |

REFCO Laboratories
10/17/96

Folder: pripl198
Protocol: priorplt

Page 26

| Line | Conc | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|---|-------|-------|--------|-------|-------|-------|-------|---|
| Sample ID: 96T14939 10-17-96 Seq: 564 11:10:18 17 Oct 1996 ICP | | | | | | | | |
| | | | 1/160 | | | | | |
| 1 | -.007 | PPM | .006 | -.008 | -.008 | -.012 | .001 | |
| 2 | -.006 | PPM | .073 | .061 | -.104 | .037 | -.019 | |
| 3 | -.000 | PPM | .000 | .000 | -.001 | -.001 | -.000 | |
| 4 | -.002 | PPM | .000 | -.002 | -.002 | -.003 | -.002 | |
| 5 | .001 | PPM | .000 | .002 | .001 | .001 | .002 | |
| 6 | .001 | PPM | .004 | -.003 | .006 | .000 | -.001 | |
| 7 | .022 | PPM | .002 | .019 | .023 | .024 | .022 | |
| 8 | .022 | PPM | .027 | .016 | .055 | -.009 | .026 | |
| 9 | .012 | PPM | .001 | .011 | .011 | .011 | .013 | |
| 10 | .020 | PPM | .017 | .002 | .042 | .017 | .019 | |
| 11 | -.011 | PPM | .067 | -.087 | -.018 | .076 | -.015 | |
| 12 | .023 | PPM | .002 | .025 | .023 | .023 | .021 | |
| 13 | .034 | PPM | .072 | .009 | -.055 | .076 | .107 | |
| 14 | -.003 | PPM | .003 | -.006 | -.002 | -.001 | -.001 | |
| 15 | .690 | PPM | .011 | .674 | .698 | .692 | .697 | |
| 16 | 1.86 | PPM | .053 | 1.78 | 1.86 | 1.87 | 1.91 | |
| Zn1 | 4.47 | PPM | .064 | 4.38 | 4.48 | 4.50 | 4.52 | |
| 17 | .028 | PPM | .015 | .048 | .011 | .026 | .027 | |
| 18 | 90.9 | PPM | 1.41 | 88.8 | 91.2 | 91.3 | 92.1 | |
| 19 | .697 | PPM | .025 | .730 | .668 | .699 | .692 | |
| 20 | .208 | PPM | .025 | .227 | .192 | .233 | .182 | |
| 21 | .093 | PPM | .390 | -.425 | .067 | .225 | .505 | |
| 22 | .014 | PPM | .015 | .006 | .027 | -.002 | .026 | |

Section V
Metals ICP Calibration Curves

E.C. Laboratories
14:59:25 16 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 1

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|------|-------|-------|--------|---|---|---|---|---|
|------|-------|-------|--------|---|---|---|---|---|

* * Peak Seq: 446 14:59:25 16 Oct 1996

XAp = 215 YAp = 141 Intensity = 7402824

* * Peak Seq: 447 15:00:02 16 Oct 1996

XAp = 215 YAp = 141 Intensity = 7409364

E.C. Laboratories
15:13:10 16 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 2

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|-----------------------|-------|-------|-------------|--------|---------|--------|----|---|
| ** Standard: 1 Rep: 1 | | | | | | | | |
| As2 | .000 | PPM | | 167 | 637 | 783 | 79 | |
| B2 | .000 | PPM | Ave. Int. = | 417 | S. D. = | 346 | | |
| | | | -2209 | -585 | 135 | -3443 | | |
| | | | Ave. Int. = | -1526 | S. D. = | 1611 | | |
| Ba1 | .000 | PPM | 108023 | 98199 | 108125 | 97389 | | |
| C1 | .000 | PPM | Ave. Int. = | 102934 | S. D. = | 5945 | | |
| | | | -1811 | -8795 | -2611 | 415 | | |
| | | | Ave. Int. = | -3201 | S. D. = | 3943 | | |
| Co1 | .000 | PPM | 631 | 1491 | 1797 | -121 | | |
| | | | Ave. Int. = | 950 | S. D. = | 868 | | |
| Cr4 | .000 | PPM | 127 | -5423 | -703 | -4505 | | |
| | | | Ave. Int. = | -2626 | S. D. = | 2747 | | |
| F1 | .000 | PPM | -643 | -241 | 373 | -2065 | | |
| | | | Ave. Int. = | -644 | S. D. = | 1035 | | |
| Mn1 | .000 | PPM | 5173 | 8307 | 19075 | 11077 | | |
| | | | Ave. Int. = | 10908 | S. D. = | 5955 | | |
| S1 | .000 | PPM | -263 | 939 | 815 | 1455 | | |
| | | | Ave. Int. = | 737 | S. D. = | 722 | | |
| Tl1 | .000 | PPM | -137 | 395 | 215 | -317 | | |
| | | | Ave. Int. = | 39 | S. D. = | 324 | | |
| Ni3 | .000 | PPM | 4509 | 517 | 2295 | 1679 | | |
| | | | Ave. Int. = | 2250 | S. D. = | 1677 | | |
| S1 | .000 | PPM | 261 | 599 | 895 | 903 | | |
| | | | Ave. Int. = | 665 | S. D. = | 304 | | |
| V2 | .000 | PPM | 141557 | 130133 | 139771 | 133773 | | |
| | | | Ave. Int. = | 136309 | S. D. = | 5295 | | |
| Ca1 | .000 | PPM | 6067 | -861 | -71 | 545 | | |
| | | | Ave. Int. = | 1420 | S. D. = | 3151 | | |
| Fe2 | .000 | PPM | 14161 | 19081 | 32729 | 16641 | | |
| | | | Ave. Int. = | 20653 | S. D. = | 8297 | | |
| Zn1 | .000 | PPM | 2533 | 797 | 2781 | 3341 | | |
| | | | Ave. Int. = | 2363 | S. D. = | 1097 | | |
| Al3 | .000 | PPM | 12111 | 7345 | 12033 | 9735 | | |
| | | | Ave. Int. = | 10306 | S. D. = | 2261 | | |
| Na2 | .000 | PPM | -2831 | -3887 | -8085 | -5895 | | |
| | | | Ave. Int. = | -5175 | S. D. = | 2319 | | |
| Ca3 | .000 | PPM | 33976 | 29000 | 39898 | 32186 | | |
| | | | Ave. Int. = | 33765 | S. D. = | 4577 | | |
| Mg4 | .000 | PPM | 24925 | 23597 | 26629 | 24237 | | |
| | | | Ave. Int. = | 24847 | S. D. = | 1306 | | |
| H1 | .000 | PPM | -1419 | -1825 | -1569 | -1927 | | |
| | | | Ave. Int. = | -1685 | S. D. = | 233 | | |
| Mo1 | .000 | PPM | 1495 | 855 | 1965 | 1557 | | |
| | | | Ave. Int. = | 1468 | S. D. = | 459 | | |

E.C. Laboratories
15:14:04 16 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 3

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|-----------------------|-------|-------|-------------|--------|----------------|--------|--------|---|
| ** Standard: 1 Rep: 2 | | | | | | | | |
| As2 | .000 | PPM | | -49 | -785 | 229 | -389 | |
| | | | Ave. Int. = | | -249 S. D. = | | 438 | |
| B2 | .000 | PPM | | -4561 | -3339 | -2315 | -2145 | |
| | | | Ave. Int. = | | -3090 S. D. = | | 1113 | |
| Ba1 | .000 | PPM | | 109989 | 105887 | 100735 | 104013 | |
| | | | Ave. Int. = | | 105156 S. D. = | | 3862 | |
| Cd1 | .000 | PPM | | -5081 | -7227 | -9067 | -8571 | |
| | | | Ave. Int. = | | -7487 S. D. = | | 1782 | |
| Co1 | .000 | PPM | | -355 | 373 | 811 | -665 | |
| | | | Ave. Int. = | | 41 S. D. = | | 673 | |
| Cr4 | .000 | PPM | | -4815 | 943 | -4537 | -10881 | |
| | | | Ave. Int. = | | -4823 S. D. = | | 4831 | |
| Fe1 | .000 | PPM | | -225 | -411 | -35 | -2219 | |
| | | | Ave. Int. = | | -723 S. D. = | | 1009 | |
| Mn1 | .000 | PPM | | 14691 | 15125 | 13517 | 10955 | |
| | | | Ave. Int. = | | 13572 S. D. = | | 1872 | |
| Mo1 | .000 | PPM | | 1179 | 313 | 713 | -173 | |
| | | | Ave. Int. = | | 508 S. D. = | | 576 | |
| Tl1 | .000 | PPM | | 169 | -425 | 31 | -385 | |
| | | | Ave. Int. = | | -153 S. D. = | | 297 | |
| Ni3 | .000 | PPM | | 5085 | 2877 | 2109 | 3051 | |
| | | | Ave. Int. = | | 3281 S. D. = | | 1271 | |
| Se1 | .000 | PPM | | 687 | 653 | 679 | 455 | |
| | | | Ave. Int. = | | 619 S. D. = | | 110 | |
| V2 | .000 | PPM | | 131645 | 133211 | 137457 | 134691 | |
| | | | Ave. Int. = | | 134251 S. D. = | | 2473 | |
| Al1 | .000 | PPM | | -867 | 373 | 1921 | 2299 | |
| | | | Ave. Int. = | | 932 S. D. = | | 1460 | |
| Fe2 | .000 | PPM | | 16521 | 22377 | 18471 | 18831 | |
| | | | Ave. Int. = | | 19050 S. D. = | | 2439 | |
| Zn1 | .000 | PPM | | 2791 | 2245 | 4327 | 2735 | |
| | | | Ave. Int. = | | 3025 S. D. = | | 902 | |
| Li3 | .000 | PPM | | 12761 | 15871 | 12049 | 13057 | |
| | | | Ave. Int. = | | 13435 S. D. = | | 1679 | |
| Na2 | .000 | PPM | | -5639 | -5221 | -3151 | -2477 | |
| | | | Ave. Int. = | | -4122 S. D. = | | 1545 | |
| Ka3 | .000 | PPM | | 37432 | 31370 | 35426 | 33218 | |
| | | | Ave. Int. = | | 34362 S. D. = | | 2634 | |
| Mg4 | .000 | PPM | | 23067 | 25261 | 24277 | 21877 | |
| | | | Ave. Int. = | | 23621 S. D. = | | 1468 | |
| Ca1 | .000 | PPM | | -727 | -1045 | -1993 | -804 | |
| | | | Ave. Int. = | | -1142 S. D. = | | 583 | |
| Mo1 | .000 | PPM | | -977 | 2411 | 1229 | 1311 | |
| | | | Ave. Int. = | | 994 S. D. = | | 1420 | |

E.C. Laboratories
15:01 16 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 4

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|------------------------|-------|-------|-------------|--------|----------------|--------|--------|---|
| *** Standard: 1 Rep: 3 | | | | | | | | |
| As2 | .000 | PPM | | 423 | 31 | -57 | -353 | |
| | | | Ave. Int. = | | 11 S. D. = | 320 | | |
| B2 | .000 | PPM | | -2491 | -2715 | -2163 | -4411 | |
| | | | Ave. Int. = | | -2945 S. D. = | 1003 | | |
| Ba1 | .000 | PPM | | 112375 | 109335 | 109853 | 98597 | |
| | | | Ave. Int. = | | 107540 S. D. = | 6108 | | |
| Cd1 | .000 | PPM | | -5611 | -6889 | -4849 | -5623 | |
| | | | Ave. Int. = | | -5743 S. D. = | 845 | | |
| Co1 | .000 | PPM | | -1499 | 125 | -65 | -203 | |
| | | | Ave. Int. = | | -411 S. D. = | 738 | | |
| Cr4 | .000 | PPM | | -1423 | 2577 | -9335 | -6927 | |
| | | | Ave. Int. = | | -3777 S. D. = | 5377 | | |
| Pb1 | .000 | PPM | | 77 | -947 | -1099 | -1211 | |
| | | | Ave. Int. = | | -795 S. D. = | 591 | | |
| Mn1 | .000 | PPM | | 10877 | 14101 | 8467 | 10267 | |
| | | | Ave. Int. = | | 10928 S. D. = | 2350 | | |
| Si1 | .000 | PPM | | -325 | 881 | 2407 | 445 | |
| | | | Ave. Int. = | | 852 S. D. = | 1150 | | |
| Tl1 | .000 | PPM | | -413 | -445 | 37 | -677 | |
| | | | Ave. Int. = | | -375 S. D. = | 298 | | |
| Li3 | .000 | PPM | | 389 | 2109 | 415 | 2445 | |
| | | | Ave. Int. = | | 1340 S. D. = | 1091 | | |
| Se1 | .000 | PPM | | -1 | 475 | 821 | 111 | |
| | | | Ave. Int. = | | 352 S. D. = | 373 | | |
| V2 | .000 | PPM | | 139917 | 136347 | 144819 | 135211 | |
| | | | Ave. Int. = | | 139074 S. D. = | 4323 | | |
| Al1 | .000 | PPM | | 2983 | 2797 | 4149 | -2671 | |
| | | | Ave. Int. = | | 1815 S. D. = | 3050 | | |
| Fe2 | .000 | PPM | | 17761 | 13127 | 18393 | 19829 | |
| | | | Ave. Int. = | | 17278 S. D. = | 2899 | | |
| Sn1 | .000 | PPM | | 3789 | 1493 | 2045 | 2327 | |
| | | | Ave. Int. = | | 2414 S. D. = | 980 | | |
| Al3 | .000 | PPM | | 9607 | 11111 | 9687 | 11455 | |
| | | | Ave. Int. = | | 10465 S. D. = | 955 | | |
| Na2 | .000 | PPM | | -5405 | -5829 | -3863 | -5415 | |
| | | | Ave. Int. = | | -5128 S. D. = | 866 | | |
| Ca3 | .000 | PPM | | 35880 | 38064 | 42650 | 34632 | |
| | | | Ave. Int. = | | 37807 S. D. = | 3527 | | |
| Mg4 | .000 | PPM | | 24677 | 20901 | 24037 | 24637 | |
| | | | Ave. Int. = | | 23563 S. D. = | 1799 | | |
| Li1 | .000 | PPM | | -1389 | -335 | -1101 | -1121 | |
| | | | Ave. Int. = | | -987 S. D. = | 454 | | |
| Mol | .000 | PPM | | 1735 | 1263 | 621 | 2183 | |
| | | | Ave. Int. = | | 1451 S. D. = | 669 | | |

E.C. Laboratories
15:18:10 16 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 5

| Name | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|------------------------|-------|-------|-------------|----------|----------|----------|----------|---|
| *** Standard: 2 Rep: 1 | | | | | | | | |
| As2 | 1.00 | PPM | | 14589 | 15423 | 14637 | 16405 | |
| | | | Ave. Int. = | 15264 | S. D. = | 852 | | |
| B2 | 1.00 | PPM | | 3825348 | 3783126 | 3784390 | 3781244 | |
| | | | Ave. Int. = | 3793527 | S. D. = | 21253 | | |
| Ba1 | 1.00 | PPM | | 8200223 | 8104945 | 8128737 | 8069307 | |
| | | | Ave. Int. = | 8125803 | S. D. = | 55299 | | |
| Ca1 | 1.00 | PPM | | 1750553 | 1737425 | 1746335 | 1723319 | |
| | | | Ave. Int. = | 1739408 | S. D. = | 12041 | | |
| Co1 | 1.00 | PPM | | 192440 | 191150 | 190944 | 190478 | |
| | | | Ave. Int. = | 191253 | S. D. = | 840 | | |
| Cr4 | 1.50 | PPM | | 2235184 | 2200806 | 2216392 | 2206742 | |
| | | | Ave. Int. = | 2214781 | S. D. = | 15042 | | |
| Fe1 | 1.00 | PPM | | 57274 | 55170 | 56058 | 57436 | |
| | | | Ave. Int. = | 56485 | S. D. = | 1071 | | |
| Mn1 | 1.50 | PPM | | 3745829 | 3704723 | 3713874 | 3693075 | |
| | | | Ave. Int. = | 3714375 | S. D. = | 22631 | | |
| Si1 | 1.00 | PPM | | 20249 | 20029 | 18675 | 19269 | |
| | | | Ave. Int. = | 19556 | S. D. = | 722 | | |
| Tl1 | 1.00 | PPM | | 4755 | 5105 | 4143 | 4561 | |
| | | | Ave. Int. = | 4641 | S. D. = | 401 | | |
| Na3 | 1.50 | PPM | | 371789 | 366485 | 371613 | 366613 | |
| | | | Ave. Int. = | 369125 | S. D. = | 2976 | | |
| Se1 | 1.00 | PPM | | 7237 | 6167 | 7383 | 7591 | |
| | | | Ave. Int. = | 7095 | S. D. = | 635 | | |
| V2 | 1.50 | PPM | | 2825548 | 2794404 | 2803506 | 2789618 | |
| | | | Ave. Int. = | 2803269 | S. D. = | 15931 | | |
| Ca1 | 1.50 | PPM | | 1090326 | 1077790 | 1076102 | 1070066 | |
| | | | Ave. Int. = | 1078571 | S. D. = | 8509 | | |
| Fe2 | 1.50 | PPM | | 1148768 | 1148320 | 1146395 | 1131904 | |
| | | | Ave. Int. = | 1143847 | S. D. = | 8028 | | |
| In1 | 1.50 | PPM | | 492242 | 485290 | 492618 | 483744 | |
| | | | Ave. Int. = | 488474 | S. D. = | 4615 | | |
| Al3 | 1.50 | PPM | | 147788 | 150484 | 153836 | 152874 | |
| | | | Ave. Int. = | 151246 | S. D. = | 2702 | | |
| Na2 | 10.0 | PPM | | 112682 | 112348 | 114852 | 111258 | |
| | | | Ave. Int. = | 112785 | S. D. = | 1506 | | |
| Ca3 | 10.0 | PPM | | 10956331 | 10806847 | 10840631 | 10772808 | |
| | | | Ave. Int. = | 10844154 | S. D. = | 79746 | | |
| Mg4 | 10.0 | PPM | | 532273 | 527065 | 529241 | 526086 | |
| | | | Ave. Int. = | 528666 | S. D. = | 2742 | | |
| Li1 | 10.0 | PPM | | 8304 | 8614 | 9010 | 8552 | |
| | | | Ave. Int. = | 8620 | S. D. = | 292 | | |
| Mo1 | 1.00 | PPM | | 77778 | 76964 | 79306 | 76498 | |
| | | | Ave. Int. = | 77637 | S. D. = | 1232 | | |

E.C. Laboratories
15:19:07 16 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 6

| Name | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|------------------------|-------|-------|-------------|----------|----------|----------|----------|---|
| *** Standard: 2 Rep: 2 | | | | | | | | |
| As2 | 1.00 | PPM | | 16661 | 16587 | 14567 | 15157 | |
| | | | Ave. Int. = | | 15743 | S. D. = | 1046 | |
| B2 | 1.00 | PPM | | 3763698 | 3731340 | 3782798 | 3783126 | |
| | | | Ave. Int. = | | 3765241 | S. D. = | 24357 | |
| Ba1 | 1.00 | PPM | | 8085885 | 7998841 | 8115059 | 8109059 | |
| | | | Ave. Int. = | | 8077211 | S. D. = | 53740 | |
| C1 | 1.00 | PPM | | 1745369 | 1724255 | 1741639 | 1744689 | |
| | | | Ave. Int. = | | 1738988 | S. D. = | 9955 | |
| Co1 | 1.00 | PPM | | 190822 | 188598 | 190190 | 189542 | |
| | | | Ave. Int. = | | 189788 | S. D. = | 950 | |
| Cr4 | 1.50 | PPM | | 2219664 | 2189811 | 2220964 | 2213712 | |
| | | | Ave. Int. = | | 2211038 | S. D. = | 14499 | |
| F1 | 1.00 | PPM | | 58594 | 58594 | 57506 | 56612 | |
| | | | Ave. Int. = | | 57827 | S. D. = | 958 | |
| Mn1 | 1.50 | PPM | | 3700123 | 3646745 | 3699547 | 3697381 | |
| | | | Ave. Int. = | | 3685949 | S. D. = | 26163 | |
| S1 | 1.00 | PPM | | 19663 | 18677 | 18469 | 19379 | |
| | | | Ave. Int. = | | 19047 | S. D. = | 566 | |
| Tl1 | 1.00 | PPM | | 4569 | 4215 | 4159 | 4519 | |
| | | | Ave. Int. = | | 4366 | S. D. = | 208 | |
| N3 | 1.50 | PPM | | 365693 | 361781 | 367087 | 369925 | |
| | | | Ave. Int. = | | 366122 | S. D. = | 3387 | |
| S1 | 1.00 | PPM | | 7447 | 7501 | 7559 | 6669 | |
| | | | Ave. Int. = | | 7294 | S. D. = | 419 | |
| V2 | 1.50 | PPM | | 2786364 | 2757700 | 2796450 | 2796266 | |
| | | | Ave. Int. = | | 2784195 | S. D. = | 18281 | |
| Cl1 | 1.50 | PPM | | 1076228 | 1066064 | 1079446 | 1082464 | |
| | | | Ave. Int. = | | 1076051 | S. D. = | 7128 | |
| Fe2 | 1.50 | PPM | | 1137102 | 1124120 | 1132014 | 1139064 | |
| | | | Ave. Int. = | | 1133075 | S. D. = | 6668 | |
| Zn1 | 1.50 | PPM | | 487466 | 483650 | 491610 | 495752 | |
| | | | Ave. Int. = | | 489620 | S. D. = | 5223 | |
| Al3 | 1.50 | PPM | | 146372 | 148516 | 140668 | 149268 | |
| | | | Ave. Int. = | | 146206 | S. D. = | 3891 | |
| Na2 | 10.0 | PPM | | 110364 | 108964 | 112548 | 109764 | |
| | | | Ave. Int. = | | 110410 | S. D. = | 1536 | |
| Ca3 | 10.0 | PPM | | 10783846 | 10716575 | 10829871 | 10812479 | |
| | | | Ave. Int. = | | 10785693 | S. D. = | 49833 | |
| Mg4 | 10.0 | PPM | | 525761 | 527185 | 527361 | 526155 | |
| | | | Ave. Int. = | | 526616 | S. D. = | 779 | |
| H1 | 10.0 | PPM | | 8342 | 8220 | 9056 | 7930 | |
| | | | Ave. Int. = | | 8387 | S. D. = | 478 | |
| Mn1 | 1.00 | PPM | | 73924 | 75100 | 76092 | 75938 | |
| | | | Ave. Int. = | | 75264 | S. D. = | 994 | |

E.C. Laboratories
15:20:04 16 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 7

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|------------------------|-------|-------|-------------|----------|----------|----------------------|----------|---|
| *** Standard: 2 Rep: 3 | | | | | | | | |
| | | | | Seq: | 453 | 15:20:04 16 Oct 1996 | ICP | |
| As2 | 1.00 | PPM | | 13719 | 13845 | 14903 | 15285 | |
| | | | Ave. Int. = | 14438 | S. D. = | 775 | | |
| B2 | 1.00 | PPM | | 3715098 | 3723684 | 3768230 | 3768438 | |
| | | | Ave. Int. = | 3743863 | S. D. = | 28474 | | |
| Ba1 | 1.00 | PPM | | 7972743 | 7992655 | 8058268 | 8042228 | |
| | | | Ave. Int. = | 8016474 | S. D. = | 40372 | | |
| C1 | 1.00 | PPM | | 1699306 | 1707913 | 1719247 | 1745703 | |
| | | | Ave. Int. = | 1718042 | S. D. = | 20168 | | |
| Co1 | 1.00 | PPM | | 185424 | 188070 | 189542 | 187726 | |
| | | | Ave. Int. = | 187691 | S. D. = | 1704 | | |
| Cr4 | 1.50 | PPM | | 2169795 | 2189523 | 2201760 | 2206872 | |
| | | | Ave. Int. = | 2191988 | S. D. = | 16489 | | |
| P1 | 1.00 | PPM | | 58794 | 58226 | 55690 | 56130 | |
| | | | Ave. Int. = | 57210 | S. D. = | 1530 | | |
| Mn1 | 1.50 | PPM | | 3640313 | 3662923 | 3686291 | 3677869 | |
| | | | Ave. Int. = | 3666849 | S. D. = | 20158 | | |
| S1 | 1.00 | PPM | | 18635 | 19277 | 17327 | 18969 | |
| | | | Ave. Int. = | 18552 | S. D. = | 858 | | |
| Tl1 | 1.00 | PPM | | 4397 | 4807 | 4571 | 4219 | |
| | | | Ave. Int. = | 4499 | S. D. = | 251 | | |
| Ni3 | 1.50 | PPM | | 361797 | 368004 | 369053 | 366175 | |
| | | | Ave. Int. = | 366257 | S. D. = | 3202 | | |
| S1 | 1.00 | PPM | | 6959 | 7261 | 7327 | 7303 | |
| | | | Ave. Int. = | 7213 | S. D. = | 171 | | |
| V2 | 1.50 | PPM | | 2727882 | 2759316 | 2780954 | 2772402 | |
| | | | Ave. Int. = | 2760139 | S. D. = | 23273 | | |
| Ca1 | 1.50 | PPM | | 1050042 | 1065314 | 1065154 | 1071468 | |
| | | | Ave. Int. = | 1062995 | S. D. = | 9122 | | |
| Fe2 | 1.50 | PPM | | 1126928 | 1124580 | 1136600 | 1136248 | |
| | | | Ave. Int. = | 1131089 | S. D. = | 6236 | | |
| In1 | 1.50 | PPM | | 479081 | 483418 | 486056 | 491154 | |
| | | | Ave. Int. = | 484927 | S. D. = | 5050 | | |
| Al3 | 1.50 | PPM | | 139554 | 138500 | 144636 | 142004 | |
| | | | Ave. Int. = | 141174 | S. D. = | 2736 | | |
| Na2 | 10.0 | PPM | | 110306 | 110032 | 114292 | 111216 | |
| | | | Ave. Int. = | 111462 | S. D. = | 1954 | | |
| Ca3 | 10.0 | PPM | | 10616430 | 10680733 | 10769326 | 10796639 | |
| | | | Ave. Int. = | 10715782 | S. D. = | 82672 | | |
| Mg4 | 10.0 | PPM | | 517489 | 521619 | 528465 | 527009 | |
| | | | Ave. Int. = | 523646 | S. D. = | 5051 | | |
| H1 | 10.0 | PPM | | 7458 | 8846 | 7654 | 8244 | |
| | | | Ave. Int. = | 8051 | S. D. = | 627 | | |
| Mo1 | 1.00 | PPM | | 75404 | 75162 | 76722 | 77762 | |
| | | | Ave. Int. = | 76263 | S. D. = | 1212 | | |

E.C. Laboratories
15:21:53 16 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 8

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|------|-------|-------|--------|---|---|---|---|---|
|------|-------|-------|--------|---|---|---|---|---|

*** Peak Seq: 454 15:21:53 16 Oct 1996

XAp = 215 YAp = 141 Intensity = 7390514

E.C. Laboratories
15:23:13 16 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 9

| Name | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|------------------------|-------|-------|-------------|----------|----------|----------|----------|---|
| *** Standard: 3 Rep: 1 | | | | | | | | |
| As2 | 3.00 | PPM | | 45581 | 44485 | 44943 | 43921 | |
| | | | Ave. Int. = | 44733 | S. D. = | 703 | | |
| B2 | 3.00 | PPM | | 11004747 | 11211419 | 11223107 | 11211421 | |
| | | | Ave. Int. = | 11162674 | S. D. = | 105428 | | |
| Ba1 | 3.00 | PPM | | 23300902 | 23737400 | 23745546 | 23702175 | |
| | | | Ave. Int. = | 23621506 | S. D. = | 214563 | | |
| C1 | 3.00 | PPM | | 5062177 | 5144463 | 5142175 | 5142481 | |
| | | | Ave. Int. = | 5122824 | S. D. = | 40444 | | |
| Co1 | 3.00 | PPM | | 552171 | 563627 | 561611 | 563613 | |
| | | | Ave. Int. = | 560256 | S. D. = | 5472 | | |
| Cr4 | 4.50 | PPM | | 6434987 | 6554676 | 6579273 | 6555948 | |
| | | | Ave. Int. = | 6531221 | S. D. = | 65145 | | |
| P1 | 3.00 | PPM | | 170185 | 171049 | 171791 | 172347 | |
| | | | Ave. Int. = | 171343 | S. D. = | 937 | | |
| Mn1 | 4.50 | PPM | | 10785645 | 10967168 | 10987218 | 10981232 | |
| | | | Ave. Int. = | 10930316 | S. D. = | 96813 | | |
| S1 | 3.00 | PPM | | 54395 | 54627 | 53509 | 54949 | |
| | | | Ave. Int. = | 54370 | S. D. = | 617 | | |
| Tl1 | 3.00 | PPM | | 13183 | 12579 | 14399 | 13557 | |
| | | | Ave. Int. = | 13430 | S. D. = | 762 | | |
| Ni3 | 4.50 | PPM | | 1061335 | 1079384 | 1082048 | 1079431 | |
| | | | Ave. Int. = | 1075550 | S. D. = | 9558 | | |
| Se1 | 3.00 | PPM | | 19429 | 20759 | 19821 | 19407 | |
| | | | Ave. Int. = | 19854 | S. D. = | 633 | | |
| V2 | 4.50 | PPM | | 7833207 | 7970735 | 7997165 | 7973431 | |
| | | | Ave. Int. = | 7943635 | S. D. = | 74570 | | |
| Ca1 | 4.50 | PPM | | 3113158 | 3176866 | 3156856 | 3171820 | |
| | | | Ave. Int. = | 3154675 | S. D. = | 28953 | | |
| Fe2 | 4.50 | PPM | | 3284452 | 3338369 | 3345036 | 3346154 | |
| | | | Ave. Int. = | 3328503 | S. D. = | 29568 | | |
| Sn1 | 4.50 | PPM | | 1416907 | 1455037 | 1451469 | 1449435 | |
| | | | Ave. Int. = | 1443212 | S. D. = | 17689 | | |
| Al3 | 4.50 | PPM | | 408202 | 416400 | 417042 | 417600 | |
| | | | Ave. Int. = | 414811 | S. D. = | 4433 | | |
| Na2 | 30.0 | PPM | | 334045 | 343981 | 343501 | 337407 | |
| | | | Ave. Int. = | 339734 | S. D. = | 4831 | | |
| Ca3 | 30.0 | PPM | | 31531656 | 32096853 | 32133220 | 32116611 | |
| | | | Ave. Int. = | 31969585 | S. D. = | 292331 | | |
| Mg4 | 30.0 | PPM | | 1509581 | 1529039 | 1526063 | 1528157 | |
| | | | Ave. Int. = | 1523210 | S. D. = | 9171 | | |
| Li1 | 30.0 | PPM | | 26362 | 26730 | 26406 | 25988 | |
| | | | Ave. Int. = | 26372 | S. D. = | 304 | | |
| Mn1 | 3.00 | PPM | | 214188 | 224754 | 225764 | 223934 | |
| | | | Ave. Int. = | 222160 | S. D. = | 5367 | | |

E.C. Laboratories
15:24:07 16 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 10

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|------------------------|-------|-------|-------------|----------|----------|----------|----------|---|
| *** Standard: 3 Rep: 2 | | | | | | | | |
| As2 | 3.00 | PPM | | 45677 | 47247 | 45197 | 43623 | |
| | | | Ave. Int. = | 45436 | S. D. = | 1492 | | |
| Fe2 | 3.00 | PPM | | 11085088 | 11130392 | 11122446 | 11070292 | |
| | | | Ave. Int. = | 11102055 | S. D. = | 28957 | | |
| Ba1 | 3.00 | PPM | | 23452716 | 23548817 | 23512118 | 23425460 | |
| | | | Ave. Int. = | 23484778 | S. D. = | 55961 | | |
| Cd1 | 3.00 | PPM | | 5107321 | 5120937 | 5113255 | 5093843 | |
| | | | Ave. Int. = | 5108839 | S. D. = | 11446 | | |
| Co1 | 3.00 | PPM | | 558997 | 556779 | 553387 | 554877 | |
| | | | Ave. Int. = | 556010 | S. D. = | 2427 | | |
| Cr4 | 4.50 | PPM | | 6486935 | 6516188 | 6535844 | 6475453 | |
| | | | Ave. Int. = | 6503605 | S. D. = | 27496 | | |
| Co1 | 3.00 | PPM | | 167945 | 170623 | 170257 | 170287 | |
| | | | Ave. Int. = | 169778 | S. D. = | 1233 | | |
| Mn1 | 4.50 | PPM | | 10874084 | 10921848 | 10888454 | 10848542 | |
| | | | Ave. Int. = | 10883232 | S. D. = | 30581 | | |
| Co1 | 3.00 | PPM | | 54255 | 54889 | 53999 | 54943 | |
| | | | Ave. Int. = | 54522 | S. D. = | 468 | | |
| Tl1 | 3.00 | PPM | | 13041 | 13851 | 13683 | 13895 | |
| | | | Ave. Int. = | 13618 | S. D. = | 395 | | |
| Li3 | 4.50 | PPM | | 1074873 | 1075091 | 1075299 | 1067032 | |
| | | | Ave. Int. = | 1073074 | S. D. = | 4032 | | |
| Se1 | 3.00 | PPM | | 21847 | 20069 | 20887 | 19463 | |
| | | | Ave. Int. = | 20567 | S. D. = | 1034 | | |
| V2 | 4.50 | PPM | | 7905218 | 7952790 | 7943796 | 7907818 | |
| | | | Ave. Int. = | 7927406 | S. D. = | 24420 | | |
| U1 | 4.50 | PPM | | 3144324 | 3156292 | 3144320 | 3122156 | |
| | | | Ave. Int. = | 3141773 | S. D. = | 14243 | | |
| Fe2 | 4.50 | PPM | | 3306542 | 3317522 | 3316124 | 3298702 | |
| | | | Ave. Int. = | 3309723 | S. D. = | 8820 | | |
| Sn1 | 4.50 | PPM | | 1431933 | 1444963 | 1444355 | 1433949 | |
| | | | Ave. Int. = | 1438800 | S. D. = | 6820 | | |
| Al3 | 4.50 | PPM | | 411776 | 416146 | 413050 | 408544 | |
| | | | Ave. Int. = | 412379 | S. D. = | 3147 | | |
| Na2 | 30.0 | PPM | | 337551 | 341135 | 335039 | 338717 | |
| | | | Ave. Int. = | 338111 | S. D. = | 2534 | | |
| Ka3 | 30.0 | PPM | | 31889924 | 31956509 | 31899508 | 31747633 | |
| | | | Ave. Int. = | 31873394 | S. D. = | 88843 | | |
| Mg4 | 30.0 | PPM | | 1514719 | 1520343 | 1516685 | 1510431 | |
| | | | Ave. Int. = | 1515545 | S. D. = | 4129 | | |
| Li1 | 30.0 | PPM | | 27008 | 25968 | 26856 | 25616 | |
| | | | Ave. Int. = | 26362 | S. D. = | 677 | | |
| Mol | 3.00 | PPM | | 220356 | 224254 | 224028 | 222830 | |
| | | | Ave. Int. = | 222867 | S. D. = | 1787 | | |

E.C. Laboratories
16:25:04 16 Oct 1996

Folder: pripl198
Protocol: priorpit

Page 11

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|------------------------|-------|-------|-------------|----------|----------|----------|----------|---|
| *** Standard: 3 Rep: 3 | | | | | | | | |
| As2 | 3.00 | PPM | | 45461 | 45373 | 45813 | 43727 | |
| | | | Ave. Int. = | 45094 | S. D. = | 931 | | |
| B2 | 3.00 | PPM | | 11177328 | 11062758 | 11016635 | 11051804 | |
| | | | Ave. Int. = | 11077131 | S. D. = | 69635 | | |
| Ba1 | 3.00 | PPM | | 23610969 | 23310385 | 23236336 | 23277142 | |
| | | | Ave. Int. = | 23358708 | S. D. = | 170879 | | |
| C11 | 3.00 | PPM | | 5149060 | 5092353 | 5083995 | 5108361 | |
| | | | Ave. Int. = | 5108442 | S. D. = | 28904 | | |
| Co1 | 3.00 | PPM | | 559533 | 553451 | 557773 | 557451 | |
| | | | Ave. Int. = | 557052 | S. D. = | 2569 | | |
| Cr4 | 4.50 | PPM | | 6566457 | 6470103 | 6461085 | 6489567 | |
| | | | Ave. Int. = | 6496803 | S. D. = | 47933 | | |
| Fe1 | 3.00 | PPM | | 170689 | 168081 | 172747 | 171667 | |
| | | | Ave. Int. = | 170796 | S. D. = | 1996 | | |
| Mn1 | 4.50 | PPM | | 10963150 | 10820387 | 10807773 | 10847292 | |
| | | | Ave. Int. = | 10859651 | S. D. = | 70941 | | |
| Se1 | 3.00 | PPM | | 56571 | 53485 | 53205 | 52421 | |
| | | | Ave. Int. = | 53921 | S. D. = | 1823 | | |
| Tl1 | 3.00 | PPM | | 13615 | 13631 | 13667 | 12707 | |
| | | | Ave. Int. = | 13405 | S. D. = | 466 | | |
| Al3 | 4.50 | PPM | | 1077303 | 1069179 | 1065907 | 1064337 | |
| | | | Ave. Int. = | 1069182 | S. D. = | 5778 | | |
| Se1 | 3.00 | PPM | | 20909 | 19519 | 20197 | 19717 | |
| | | | Ave. Int. = | 20086 | S. D. = | 618 | | |
| V2 | 4.50 | PPM | | 7988461 | 7883335 | 7874265 | 7884297 | |
| | | | Ave. Int. = | 7907590 | S. D. = | 54103 | | |
| Ca1 | 4.50 | PPM | | 3159830 | 3120690 | 3116948 | 3123502 | |
| | | | Ave. Int. = | 3130243 | S. D. = | 19907 | | |
| Fe2 | 4.50 | PPM | | 3343186 | 3302996 | 3300756 | 3297014 | |
| | | | Ave. Int. = | 3310988 | S. D. = | 21607 | | |
| Sn1 | 4.50 | PPM | | 1448821 | 1439557 | 1428899 | 1430979 | |
| | | | Ave. Int. = | 1437064 | S. D. = | 9095 | | |
| Al3 | 4.50 | PPM | | 420338 | 411082 | 419536 | 409850 | |
| | | | Ave. Int. = | 415202 | S. D. = | 5501 | | |
| Na2 | 30.0 | PPM | | 337455 | 340087 | 339223 | 338325 | |
| | | | Ave. Int. = | 338773 | S. D. = | 1135 | | |
| As3 | 30.0 | PPM | | 32117982 | 31765855 | 31638971 | 31744345 | |
| | | | Ave. Int. = | 31816788 | S. D. = | 208310 | | |
| Mg4 | 30.0 | PPM | | 1535543 | 1513277 | 1507607 | 1516303 | |
| | | | Ave. Int. = | 1518183 | S. D. = | 12122 | | |
| Li1 | 30.0 | PPM | | 26858 | 26556 | 26552 | 26714 | |
| | | | Ave. Int. = | 26670 | S. D. = | 146 | | |
| Mn1 | 3.00 | PPM | | 224836 | 219644 | 221894 | 222620 | |
| | | | Ave. Int. = | 222249 | S. D. = | 2140 | | |

E.C. Laboratories
15:28:12 16 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 12

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|------------------------|-------|-------|-------------|----------|----------|----------|----------|---|
| *** Standard: 4 Rep: 1 | | | | | | | | |
| As2 | 7.00 | PPM | | 108797 | 106352 | 106522 | 105938 | |
| | | | Ave. Int. = | 106902 | S. D. = | 1287 | | |
| Br2 | 7.00 | PPM | | 26587756 | 26539356 | 26474251 | 26325401 | |
| | | | Ave. Int. = | 26481691 | S. D. = | 114101 | | |
| Ba1 | 7.00 | PPM | | 56084857 | 56021756 | 55970140 | 55525157 | |
| | | | Ave. Int. = | 55900478 | S. D. = | 254573 | | |
| Cl1 | 7.00 | PPM | | 11984908 | 11951291 | 11905750 | 11859182 | |
| | | | Ave. Int. = | 11925283 | S. D. = | 54719 | | |
| Co1 | 7.00 | PPM | | 1315501 | 1311885 | 1311087 | 1303855 | |
| | | | Ave. Int. = | 1310582 | S. D. = | 4879 | | |
| Cr4 | 10.5 | PPM | | 15428280 | 15404388 | 15375053 | 15265523 | |
| | | | Ave. Int. = | 15368311 | S. D. = | 71900 | | |
| Fe1 | 7.00 | PPM | | 403047 | 401071 | 397957 | 399309 | |
| | | | Ave. Int. = | 400346 | S. D. = | 2206 | | |
| Mn1 | 10.5 | PPM | | 25982435 | 25901780 | 25815934 | 25699311 | |
| | | | Ave. Int. = | 25849865 | S. D. = | 121227 | | |
| Si1 | 7.00 | PPM | | 128319 | 129483 | 128845 | 129691 | |
| | | | Ave. Int. = | 129085 | S. D. = | 624 | | |
| Tl1 | 7.00 | PPM | | 31831 | 32569 | 31559 | 32703 | |
| | | | Ave. Int. = | 32166 | S. D. = | 557 | | |
| Ni3 | 10.5 | PPM | | 2514405 | 2509757 | 2503373 | 2495093 | |
| | | | Ave. Int. = | 2505657 | S. D. = | 8370 | | |
| Se1 | 7.00 | PPM | | 46839 | 48813 | 47463 | 46781 | |
| | | | Ave. Int. = | 47474 | S. D. = | 945 | | |
| V2 | 10.5 | PPM | | 18715540 | 18682637 | 18580335 | 18495715 | |
| | | | Ave. Int. = | 18618557 | S. D. = | 100105 | | |
| Ca1 | 10.5 | PPM | | 7455551 | 7449737 | 7429558 | 7393680 | |
| | | | Ave. Int. = | 7432132 | S. D. = | 27950 | | |
| Fe2 | 10.5 | PPM | | 7854768 | 7824686 | 7810438 | 7779680 | |
| | | | Ave. Int. = | 7817393 | S. D. = | 31202 | | |
| Zn1 | 10.5 | PPM | | 3408584 | 3400184 | 3385935 | 3380405 | |
| | | | Ave. Int. = | 3393777 | S. D. = | 12918 | | |
| Al3 | 10.5 | PPM | | 983082 | 979232 | 972322 | 971026 | |
| | | | Ave. Int. = | 976416 | S. D. = | 5721 | | |
| Na2 | 70.0 | PPM | | 822242 | 817762 | 822034 | 814394 | |
| | | | Ave. Int. = | 819108 | S. D. = | 3760 | | |
| Ca3 | 70.0 | PPM | | 76013041 | 75759075 | 75428004 | 75114858 | |
| | | | Ave. Int. = | 75578745 | S. D. = | 391172 | | |
| Mg4 | 70.0 | PPM | | 3615192 | 3604448 | 3579594 | 3564700 | |
| | | | Ave. Int. = | 3590984 | S. D. = | 23006 | | |
| Li1 | 70.0 | PPM | | 64695 | 65433 | 65591 | 64639 | |
| | | | Ave. Int. = | 65090 | S. D. = | 493 | | |
| Mo1 | 7.00 | PPM | | 533665 | 524971 | 528811 | 523915 | |
| | | | Ave. Int. = | 527841 | S. D. = | 4416 | | |

E.C. Laboratories
15:29:09 16 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 13

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|------------------------|-------|-------|-------------|----------|----------|----------------------|----------|---|
| *** Standard: 4 Rep: 2 | | | | | | | | |
| | | | | Seq: | 459 | 15:29:09 16 Oct 1996 | ICP | |
| As2 | 7.00 | PPM | | 104888 | 107586 | 102370 | 105456 | |
| | | | Ave. Int. = | 105075 | S. D. = | 2145 | | |
| B2 | 7.00 | PPM | | 26128356 | 26179169 | 26274193 | 26339183 | |
| | | | Ave. Int. = | 26230225 | S. D. = | 94497 | | |
| Ba1 | 7.00 | PPM | | 55216768 | 55276650 | 55485743 | 55613296 | |
| | | | Ave. Int. = | 55398114 | S. D. = | 184049 | | |
| C1 | 7.00 | PPM | | 11777905 | 11810381 | 11869612 | 11878700 | |
| | | | Ave. Int. = | 11834150 | S. D. = | 48203 | | |
| Co1 | 7.00 | PPM | | 1295903 | 1298421 | 1303541 | 1306843 | |
| | | | Ave. Int. = | 1301177 | S. D. = | 4936 | | |
| Cr4 | 10.5 | PPM | | 15166721 | 15203944 | 15229964 | 15310659 | |
| | | | Ave. Int. = | 15227822 | S. D. = | 61019 | | |
| P1 | 7.00 | PPM | | 398933 | 399015 | 398767 | 400991 | |
| | | | Ave. Int. = | 399427 | S. D. = | 1048 | | |
| Mn1 | 10.5 | PPM | | 25619570 | 25598090 | 25704328 | 25727938 | |
| | | | Ave. Int. = | 25662482 | S. D. = | 63307 | | |
| S1 | 7.00 | PPM | | 129393 | 126343 | 128823 | 130403 | |
| | | | Ave. Int. = | 128741 | S. D. = | 1727 | | |
| Tl1 | 7.00 | PPM | | 31887 | 31043 | 32575 | 31279 | |
| | | | Ave. Int. = | 31696 | S. D. = | 685 | | |
| Ni3 | 10.5 | PPM | | 2479003 | 2480973 | 2503645 | 2496277 | |
| | | | Ave. Int. = | 2489975 | S. D. = | 11944 | | |
| S1 | 7.00 | PPM | | 46973 | 47797 | 47284 | 48661 | |
| | | | Ave. Int. = | 47679 | S. D. = | 738 | | |
| V2 | 10.5 | PPM | | 18385019 | 18410920 | 18486545 | 18545373 | |
| | | | Ave. Int. = | 18456964 | S. D. = | 73001 | | |
| Cl1 | 10.5 | PPM | | 7348111 | 7347845 | 7368871 | 7401304 | |
| | | | Ave. Int. = | 7366533 | S. D. = | 25187 | | |
| Fe2 | 10.5 | PPM | | 7757184 | 7750614 | 7772248 | 7781160 | |
| | | | Ave. Int. = | 7765302 | S. D. = | 13921 | | |
| Zn1 | 10.5 | PPM | | 3358525 | 3373937 | 3388589 | 3383351 | |
| | | | Ave. Int. = | 3376101 | S. D. = | 13192 | | |
| Al3 | 10.5 | PPM | | 963962 | 969824 | 970106 | 977466 | |
| | | | Ave. Int. = | 970340 | S. D. = | 5531 | | |
| Na2 | 70.0 | PPM | | 807560 | 817568 | 815560 | 816512 | |
| | | | Ave. Int. = | 814300 | S. D. = | 4568 | | |
| Ca3 | 70.0 | PPM | | 74568695 | 74715853 | 75048276 | 75101592 | |
| | | | Ave. Int. = | 74858604 | S. D. = | 257840 | | |
| Mg4 | 70.0 | PPM | | 3547112 | 3550420 | 3567218 | 3575706 | |
| | | | Ave. Int. = | 3560114 | S. D. = | 13621 | | |
| H1 | 70.0 | PPM | | 63582 | 64821 | 66185 | 64869 | |
| | | | Ave. Int. = | 64864 | S. D. = | 1063 | | |
| Mo1 | 7.00 | PPM | | 522715 | 521777 | 524107 | 528649 | |
| | | | Ave. Int. = | 524312 | S. D. = | 3046 | | |

E.C. Laboratories
15:30:06 16 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 14

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|------------------------|-------|-------|-------------|----------|----------|----------|----------|---|
| *** Standard: 4 Rep: 3 | | | | | | | | |
| As2 | 7.00 | PPM | | 108653 | 107115 | 107739 | 104386 | |
| | | | Ave. Int. = | 106973 | S. D. = | 1837 | | |
| Ba2 | 7.00 | PPM | | 26271084 | 26097164 | 26461315 | 26150791 | |
| | | | Ave. Int. = | 26245089 | S. D. = | 161455 | | |
| Ba1 | 7.00 | PPM | | 55414854 | 55044011 | 55730009 | 55089005 | |
| | | | Ave. Int. = | 55319470 | S. D. = | 319704 | | |
| Ca1 | 7.00 | PPM | | 11866660 | 11784905 | 11974211 | 11818423 | |
| | | | Ave. Int. = | 11861050 | S. D. = | 82567 | | |
| Co1 | 7.00 | PPM | | 1303637 | 1294551 | 1313815 | 1300324 | |
| | | | Ave. Int. = | 1303082 | S. D. = | 8081 | | |
| Cr4 | 10.5 | PPM | | 15241210 | 15165377 | 15371227 | 15206368 | |
| | | | Ave. Int. = | 15246046 | S. D. = | 89023 | | |
| Fe1 | 7.00 | PPM | | 401407 | 395503 | 400423 | 396967 | |
| | | | Ave. Int. = | 398575 | S. D. = | 2796 | | |
| Mn1 | 10.5 | PPM | | 25705823 | 25536658 | 25836242 | 25537690 | |
| | | | Ave. Int. = | 25654103 | S. D. = | 145138 | | |
| Si1 | 7.00 | PPM | | 128407 | 125845 | 130301 | 128953 | |
| | | | Ave. Int. = | 128377 | S. D. = | 1866 | | |
| Tl1 | 7.00 | PPM | | 31137 | 32989 | 31574 | 32025 | |
| | | | Ave. Int. = | 31931 | S. D. = | 793 | | |
| Ni3 | 10.5 | PPM | | 2497491 | 2474461 | 2505813 | 2481421 | |
| | | | Ave. Int. = | 2489797 | S. D. = | 14388 | | |
| Si1 | 7.00 | PPM | | 47021 | 45999 | 48931 | 48009 | |
| | | | Ave. Int. = | 47490 | S. D. = | 1263 | | |
| V2 | 10.5 | PPM | | 18511737 | 18383499 | 18623985 | 18411246 | |
| | | | Ave. Int. = | 18482617 | S. D. = | 109165 | | |
| Cr1 | 10.5 | PPM | | 7382644 | 7325817 | 7424368 | 7329671 | |
| | | | Ave. Int. = | 7365625 | S. D. = | 46967 | | |
| Fe2 | 10.5 | PPM | | 7792830 | 7732846 | 7844392 | 7764174 | |
| | | | Ave. Int. = | 7783561 | S. D. = | 47379 | | |
| Zn1 | 10.5 | PPM | | 3388231 | 3363060 | 3397182 | 3347477 | |
| | | | Ave. Int. = | 3373988 | S. D. = | 22826 | | |
| Al3 | 10.5 | PPM | | 971232 | 960464 | 979312 | 966000 | |
| | | | Ave. Int. = | 969252 | S. D. = | 8019 | | |
| Na2 | 70.0 | PPM | | 816600 | 810128 | 819426 | 809546 | |
| | | | Ave. Int. = | 813925 | S. D. = | 4865 | | |
| Cl3 | 70.0 | PPM | | 74993786 | 74478678 | 75641235 | 74588821 | |
| | | | Ave. Int. = | 74925630 | S. D. = | 525973 | | |
| Mg4 | 70.0 | PPM | | 3565978 | 3546832 | 3591798 | 3549324 | |
| | | | Ave. Int. = | 3563483 | S. D. = | 20702 | | |
| H1 | 70.0 | PPM | | 64615 | 63196 | 66869 | 64219 | |
| | | | Ave. Int. = | 64725 | S. D. = | 1549 | | |
| Mo1 | 7.00 | PPM | | 526081 | 521121 | 531113 | 524451 | |
| | | | Ave. Int. = | 525692 | S. D. = | 4162 | | |

E.C. Laboratories
1:33:15 16 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 15

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|------------------------|-------|-------|-------------|-----------|-----------|-----------|-----------|---|
| *** Standard: 5 Rep: 1 | | | | | | | | |
| As2 | 10.0 | PPM | | 148795 | 150307 | 148131 | 151501 | |
| | | | Ave. Int. = | 149684 | S. D. = | 1516 | | |
| B2 | 10.0 | PPM | | 36522638 | 36766913 | 36879091 | 36894907 | |
| | | | Ave. Int. = | 36765887 | S. D. = | 171884 | | |
| Ba1 | 10.0 | PPM | | 76700200 | 77192513 | 77285102 | 77450370 | |
| | | | Ave. Int. = | 77157046 | S. D. = | 322699 | | |
| Ca1 | 10.0 | PPM | | 16392276 | 16528760 | 16558811 | 16583947 | |
| | | | Ave. Int. = | 16515949 | S. D. = | 85479 | | |
| Co1 | 10.0 | PPM | | 1804555 | 1817081 | 1822875 | 1822281 | |
| | | | Ave. Int. = | 1816698 | S. D. = | 8503 | | |
| Cr4 | 15.0 | PPM | | 21207269 | 21307630 | 21329158 | 21381366 | |
| | | | Ave. Int. = | 21306356 | S. D. = | 72953 | | |
| Fe1 | 10.0 | PPM | | 554167 | 557509 | 560541 | 557775 | |
| | | | Ave. Int. = | 557498 | S. D. = | 2610 | | |
| Mn1 | 15.0 | PPM | | 35659964 | 35939455 | 35947044 | 35989614 | |
| | | | Ave. Int. = | 35884019 | S. D. = | 150993 | | |
| Se1 | 10.0 | PPM | | 180454 | 180788 | 178988 | 181704 | |
| | | | Ave. Int. = | 180484 | S. D. = | 1128 | | |
| Tl1 | 10.0 | PPM | | 44513 | 45433 | 43973 | 45785 | |
| | | | Ave. Int. = | 44926 | S. D. = | 831 | | |
| Ni3 | 15.0 | PPM | | 3455852 | 3481087 | 3473053 | 3498127 | |
| | | | Ave. Int. = | 3477030 | S. D. = | 17568 | | |
| Si1 | 10.0 | PPM | | 66122 | 66788 | 66890 | 67546 | |
| | | | Ave. Int. = | 66837 | S. D. = | 583 | | |
| V2 | 15.0 | PPM | | 25651199 | 25832532 | 25858045 | 25908461 | |
| | | | Ave. Int. = | 25812559 | S. D. = | 112104 | | |
| Cr1 | 15.0 | PPM | | 10215881 | 10292426 | 10305909 | 10336286 | |
| | | | Ave. Int. = | 10287626 | S. D. = | 51226 | | |
| Fe2 | 15.0 | PPM | | 10815527 | 10908027 | 10922051 | 10935037 | |
| | | | Ave. Int. = | 10895161 | S. D. = | 54223 | | |
| Sn1 | 15.0 | PPM | | 4691061 | 4742132 | 4748052 | 4740710 | |
| | | | Ave. Int. = | 4730489 | S. D. = | 26477 | | |
| Al3 | 15.0 | PPM | | 1349814 | 1358280 | 1366694 | 1377587 | |
| | | | Ave. Int. = | 1363094 | S. D. = | 11868 | | |
| Na2 | 100. | PPM | | 1136774 | 1146950 | 1142892 | 1149332 | |
| | | | Ave. Int. = | 1143987 | S. D. = | 5495 | | |
| Ca3 | 100. | PPM | | 104376505 | 105052425 | 105087756 | 105314332 | |
| | | | Ave. Int. = | 104957755 | S. D. = | 404500 | | |
| Mg4 | 100. | PPM | | 5009806 | 5035478 | 5049968 | 5052191 | |
| | | | Ave. Int. = | 5036861 | S. D. = | 19499 | | |
| Li1 | 100. | PPM | | 92505 | 90377 | 92021 | 91791 | |
| | | | Ave. Int. = | 91674 | S. D. = | 914 | | |
| Mn1 | 10.0 | PPM | | 733961 | 738194 | 738496 | 736938 | |
| | | | Ave. Int. = | 736897 | S. D. = | 2070 | | |

E.C. Laboratories
15:34:12 16 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 16

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|------------------------|-------|-------|--------|-------------|-----------|-----------|-----------|---|
| *** Standard: 5 Rep: 2 | | | | | | | | |
| A52 | 10.0 | PPM | | 150191 | 148937 | 148077 | 152029 | |
| | | | | Ave. Int. = | 149809 | S. D. = | 1716 | |
| B52 | 10.0 | PPM | | 36591460 | 36462554 | 36611540 | 37195395 | |
| | | | | Ave. Int. = | 36715237 | S. D. = | 326841 | |
| Ba1 | 10.0 | PPM | | 76759484 | 76444000 | 76674920 | 78064166 | |
| | | | | Ave. Int. = | 76985643 | S. D. = | 731274 | |
| C1 | 10.0 | PPM | | 16482422 | 16406100 | 16448982 | 16698511 | |
| | | | | Ave. Int. = | 16509004 | S. D. = | 130143 | |
| Co1 | 10.0 | PPM | | 1811707 | 1807265 | 1811249 | 1839881 | |
| | | | | Ave. Int. = | 1817526 | S. D. = | 15037 | |
| C4 | 15.0 | PPM | | 21219701 | 21125944 | 21225650 | 21607032 | |
| | | | | Ave. Int. = | 21294582 | S. D. = | 213247 | |
| P1 | 10.0 | PPM | | 554303 | 558871 | 557493 | 563093 | |
| | | | | Ave. Int. = | 558440 | S. D. = | 3645 | |
| Mn1 | 15.0 | PPM | | 35774510 | 35605869 | 35781222 | 36365031 | |
| | | | | Ave. Int. = | 35881658 | S. D. = | 332304 | |
| S1 | 10.0 | PPM | | 179900 | 179424 | 179752 | 181732 | |
| | | | | Ave. Int. = | 180202 | S. D. = | 1039 | |
| Tl1 | 10.0 | PPM | | 45183 | 45523 | 44843 | 45859 | |
| | | | | Ave. Int. = | 45352 | S. D. = | 437 | |
| N3 | 15.0 | PPM | | 3470277 | 3442682 | 3466779 | 3526475 | |
| | | | | Ave. Int. = | 3476553 | S. D. = | 35470 | |
| S1 | 10.0 | PPM | | 65610 | 67538 | 65776 | 67522 | |
| | | | | Ave. Int. = | 66612 | S. D. = | 1063 | |
| V2 | 15.0 | PPM | | 25667311 | 25602400 | 25719855 | 26133060 | |
| | | | | Ave. Int. = | 25780657 | S. D. = | 239797 | |
| C1 | 15.0 | PPM | | 10252082 | 10204291 | 10240730 | 10410984 | |
| | | | | Ave. Int. = | 10277022 | S. D. = | 91606 | |
| Fe2 | 15.0 | PPM | | 10841855 | 10796391 | 10857650 | 11037116 | |
| | | | | Ave. Int. = | 10883253 | S. D. = | 105811 | |
| Z1 | 15.0 | PPM | | 4715268 | 4698500 | 4711374 | 4791100 | |
| | | | | Ave. Int. = | 4729061 | S. D. = | 41976 | |
| Al3 | 15.0 | PPM | | 1356256 | 1351502 | 1354654 | 1380677 | |
| | | | | Ave. Int. = | 1360772 | S. D. = | 13416 | |
| Na2 | 100. | PPM | | 1139424 | 1137330 | 1138763 | 1152492 | |
| | | | | Ave. Int. = | 1142002 | S. D. = | 7048 | |
| O3 | 100. | PPM | | 104522138 | 103985248 | 104492945 | 106188339 | |
| | | | | Ave. Int. = | 104797168 | S. D. = | 959646 | |
| Mg4 | 100. | PPM | | 5027832 | 4997235 | 5028486 | 5106316 | |
| | | | | Ave. Int. = | 5039967 | S. D. = | 46574 | |
| H1 | 100. | PPM | | 90935 | 90579 | 91371 | 92213 | |
| | | | | Ave. Int. = | 91275 | S. D. = | 705 | |
| Mo1 | 10.0 | PPM | | 733715 | 731267 | 736075 | 747088 | |
| | | | | Ave. Int. = | 737036 | S. D. = | 6983 | |

E.C. Laboratories
15:35:09 16 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 17

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|------------------------|-------|-------|-------------|-----------|-----------|-----------|-----------|---|
| *** Standard: 5 Rep: 3 | | | | | | | | |
| As2 | 10.0 | PPM | | 149741 | 152973 | 151749 | 149157 | |
| | | | Ave. Int. = | 150905 | S. D. = | 1770 | | |
| B2 | 10.0 | PPM | | 36761655 | 36956561 | 36769585 | 36383366 | |
| | | | Ave. Int. = | 36717792 | S. D. = | 240456 | | |
| Ba1 | 10.0 | PPM | | 77063957 | 77470372 | 76889190 | 76246612 | |
| | | | Ave. Int. = | 76917533 | S. D. = | 509251 | | |
| C1 | 10.0 | PPM | | 16534843 | 16643598 | 16570955 | 16368330 | |
| | | | Ave. Int. = | 16529432 | S. D. = | 116535 | | |
| Co1 | 10.0 | PPM | | 1819443 | 1825801 | 1818945 | 1804155 | |
| | | | Ave. Int. = | 1817086 | S. D. = | 9168 | | |
| Cr4 | 15.0 | PPM | | 21305232 | 21466658 | 21353420 | 21120006 | |
| | | | Ave. Int. = | 21311329 | S. D. = | 144384 | | |
| Fe1 | 10.0 | PPM | | 560639 | 563973 | 563279 | 556573 | |
| | | | Ave. Int. = | 561116 | S. D. = | 3352 | | |
| Mn1 | 15.0 | PPM | | 35948488 | 36130291 | 35943254 | 35625846 | |
| | | | Ave. Int. = | 35911970 | S. D. = | 209637 | | |
| S1 | 10.0 | PPM | | 181554 | 182278 | 178522 | 180464 | |
| | | | Ave. Int. = | 180705 | S. D. = | 1635 | | |
| Tl1 | 10.0 | PPM | | 44339 | 43631 | 43777 | 43469 | |
| | | | Ave. Int. = | 43804 | S. D. = | 378 | | |
| N3 | 15.0 | PPM | | 3480693 | 3504471 | 3481965 | 3460892 | |
| | | | Ave. Int. = | 3482005 | S. D. = | 17816 | | |
| S1 | 10.0 | PPM | | 66940 | 66916 | 66522 | 66930 | |
| | | | Ave. Int. = | 66827 | S. D. = | 204 | | |
| V2 | 15.0 | PPM | | 25861885 | 25976126 | 25859005 | 25574513 | |
| | | | Ave. Int. = | 25817882 | S. D. = | 171170 | | |
| Cl1 | 15.0 | PPM | | 10297870 | 10349604 | 10268136 | 10172781 | |
| | | | Ave. Int. = | 10272098 | S. D. = | 74276 | | |
| Fe2 | 15.0 | PPM | | 10906010 | 10992944 | 10916491 | 10814885 | |
| | | | Ave. Int. = | 10907583 | S. D. = | 72941 | | |
| Zr1 | 15.0 | PPM | | 4743164 | 4766076 | 4733518 | 4703076 | |
| | | | Ave. Int. = | 4736459 | S. D. = | 26110 | | |
| Al3 | 15.0 | PPM | | 1361910 | 1371459 | 1364093 | 1353312 | |
| | | | Ave. Int. = | 1362694 | S. D. = | 7470 | | |
| Na2 | 100. | PPM | | 1149612 | 1158198 | 1144494 | 1133684 | |
| | | | Ave. Int. = | 1146497 | S. D. = | 10244 | | |
| Cl3 | 100. | PPM | | 105023691 | 105564140 | 104953866 | 103853233 | |
| | | | Ave. Int. = | 104848733 | S. D. = | 717517 | | |
| Mg4 | 100. | PPM | | 5059190 | 5063288 | 5062356 | 5005409 | |
| | | | Ave. Int. = | 5047561 | S. D. = | 28156 | | |
| H1 | 100. | PPM | | 91549 | 92237 | 90809 | 91011 | |
| | | | Ave. Int. = | 91402 | S. D. = | 639 | | |
| Mo1 | 10.0 | PPM | | 742928 | 744842 | 736507 | 734401 | |
| | | | Ave. Int. = | 739670 | S. D. = | 5004 | | |

E.C. Laboratories
15:40:48 16 Oct 1996

Folder: pripl198
Protocol: priorpl1

Page 18

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|------------------------|-----------|-------------|-------------|----------|----------|-------------|-----|---|
| *** Standard: 1 Rep: 1 | | | | Seq: 464 | 15:40:48 | 16 Oct 1996 | ICP | |
| Ag1 | .000 | PPM | 2269 | 2101 | 527 | 805 | | |
| | | | Ave. Int. = | 1426 | S. D. = | 887 | | |
| *** Standard: 1 Rep: 2 | | | | Seq: 465 | 15:41:36 | 16 Oct 1996 | ICP | |
| Ag1 | .000 | PPM | 1531 | -1011 | -5071 | -1601 | | |
| | | | Ave. Int. = | -1538 | S. D. = | 2719 | | |
| *** Standard: 1 Rep: 3 | | | | Seq: 466 | 15:42:24 | 16 Oct 1996 | ICP | |
| Ag1 | .000 | PPM | 1335 | -1231 | 2393 | 1291 | | |
| | | | Ave. Int. = | 947 | S. D. = | 1539 | | |
| ** Peak | | | | Seq: 467 | 15:44:04 | 16 Oct 1996 | | |
| XAp = 215 | YAp = 141 | Intensity = | | 7593713 | | | | |
| *** Standard: 5 Rep: 1 | | | | Seq: 468 | 15:45:24 | 16 Oct 1996 | ICP | |
| Ag1 | 1.00 | PPM | 319258 | 302288 | 285838 | 279594 | | |
| | | | Ave. Int. = | 296745 | S. D. = | 17801 | | |
| *** Standard: 5 Rep: 2 | | | | Seq: 469 | 15:46:12 | 16 Oct 1996 | ICP | |
| Ag1 | 1.00 | PPM | 269404 | 266626 | 269376 | 250385 | | |
| | | | Ave. Int. = | 263948 | S. D. = | 9135 | | |
| *** Standard: 5 Rep: 3 | | | | Seq: 470 | 15:47:00 | 16 Oct 1996 | ICP | |
| Ag1 | 1.00 | PPM | 240149 | 239381 | 238911 | 229109 | | |
| | | | Ave. Int. = | 236888 | S. D. = | 5211 | | |
| *** Standard: 6 Rep: 1 | | | | Seq: 471 | 15:50:01 | 16 Oct 1996 | ICP | |
| Ag1 | 2.00 | PPM | 590292 | 565190 | 557784 | 517367 | | |
| | | | Ave. Int. = | 557658 | S. D. = | 30249 | | |
| *** Standard: 6 Rep: 2 | | | | Seq: 472 | 15:50:49 | 16 Oct 1996 | ICP | |
| Ag1 | 2.00 | PPM | 523155 | 551901 | 540687 | 553615 | | |
| | | | Ave. Int. = | 542340 | S. D. = | 14016 | | |
| *** Standard: 6 Rep: 3 | | | | Seq: 473 | 15:51:37 | 16 Oct 1996 | ICP | |
| Ag1 | 2.00 | PPM | 559132 | 556832 | 560784 | 540109 | | |
| | | | Ave. Int. = | 554214 | S. D. = | 9542 | | |

E.C. Laboratories
15:54:38 16 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 19

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|---------------------------|-------|-------|--------|------|-------|----------|-------------|-----|
| *** Update Intercept | | | | | | | | |
| | | | Corr. | Seq: | 474 | 15:54:38 | 16 Oct 1996 | ICP |
| Ag1 | .000 | PPM | -.016 | Avg. | .016 | .010 | .016 | |
| As2 | .000 | PPM | .022 | | -.022 | .027 | -.022 | |
| Br2 | .000 | PPM | .021 | | -.021 | .000 | -.021 | |
| Br1 | .000 | PPM | .028 | | -.028 | .001 | -.028 | |
| Cd1 | .000 | PPM | .048 | | -.048 | .000 | -.048 | |
| Cl1 | .000 | PPM | .039 | | -.039 | .006 | -.039 | |
| Cl4 | .000 | PPM | .051 | | -.051 | .003 | -.051 | |
| Pb1 | .000 | PPM | .041 | | -.041 | .031 | -.041 | |
| Mn1 | .000 | PPM | .037 | | -.037 | .001 | -.037 | |
| Si1 | .000 | PPM | -.003 | | .003 | .027 | .003 | |
| Tl1 | .000 | PPM | -.003 | | .003 | .061 | .003 | |
| Ni3 | .000 | PPM | .067 | | -.067 | .004 | -.067 | |
| S11 | .000 | PPM | -.025 | | .025 | .033 | .025 | |
| V2 | .000 | PPM | .027 | | -.027 | .001 | -.027 | |
| Cu1 | .000 | PPM | .045 | | -.045 | .004 | -.045 | |
| F2 | .000 | PPM | .039 | | -.039 | .003 | -.039 | |
| Zn1 | .000 | PPM | .033 | | -.033 | .001 | -.033 | |
| Al3 | .000 | PPM | -.005 | | .005 | .047 | .005 | |
| Na2 | .000 | PPM | .021 | | -.021 | .042 | -.021 | |
| Cl3 | .000 | PPM | .241 | | -.241 | .003 | -.241 | |
| Ng4 | .000 | PPM | -.017 | | .017 | .023 | .017 | |
| K1 | .000 | PPM | -.080 | | .080 | .378 | .080 | |
| Mo1 | .000 | PPM | .005 | | -.005 | .009 | -.005 | |
| *** Check Standard: 2 Ck2 | | | | | | | | |
| Line | Flag | %Rcv. | Found | True | Units | SD/RSD | | |
| Ag1 | H | 136. | 1.36 | 1.00 | PPM | .058 | | |
| Ba1 | | 107. | 1.07 | 1.00 | PPM | .017 | | |
| Al3 | | 107. | 1.07 | 1.00 | PPM | .028 | | |

E.C. Laboratories
10:02:44 16 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 20

E.C. Laboratories
6:10:50 16 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 21

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|---------------------------|------------------|-------|---------|-------|--------|----------|-------------|-----|
| ----- | | | | | | | | |
| * | Update Intercept | | | Seq: | 479 | 16:10:50 | 16 Oct 1996 | ICP |
| | | | | Corr. | Avg. | SD/RSD | Old | |
| Ag1 | .000 | PPM | .001 | -.001 | .002 | | -.001 | |
| As2 | .000 | PPM | .005 | -.005 | .026 | | .016 | |
| Ba2 | .000 | PPM | .021 | -.021 | .000 | | .000 | |
| Ba1 | .000 | PPM | .023 | -.023 | .001 | | .005 | |
| Cd1 | .000 | PPM | .049 | -.049 | .001 | | -.001 | |
| Cr1 | .000 | PPM | .039 | -.039 | .004 | | .000 | |
| Cr4 | .000 | PPM | .053 | -.053 | .003 | | -.002 | |
| Pb1 | .000 | PPM | .062 | -.062 | .019 | | -.022 | |
| Vn1 | .000 | PPM | .038 | -.038 | .001 | | -.001 | |
| Sr1 | .000 | PPM | .009 | -.009 | .036 | | -.012 | |
| Tl1 | .000 | PPM | .033 | -.033 | .064 | | -.036 | |
| Ni3 | .000 | PPM | .072 | -.072 | .007 | | -.072 | |
| Si1 | .000 | PPM | -.028 | .028 | .041 | | .003 | |
| Vn2 | .000 | PPM | .034 | -.034 | .001 | | -.008 | |
| Cu1 | .000 | PPM | .051 | -.051 | .003 | | -.006 | |
| F12 | .000 | PPM | .033 | -.033 | .004 | | .006 | |
| Zn1 | .000 | PPM | .037 | -.037 | .004 | | -.004 | |
| Al3 | .000 | PPM | -.011 | .011 | .034 | | .007 | |
| Na2 | .000 | PPM | .052 | -.052 | .034 | | -.032 | |
| Cl3 | .000 | PPM | .234 | -.234 | .002 | | .007 | |
| Mg4 | .000 | PPM | -.044 | .044 | .024 | | .027 | |
| K1 | .000 | PPM | .173 | -.173 | .759 | | -.254 | |
| Mg1 | .000 | PPM | -.010 | .010 | .016 | | .015 | |
| ----- | | | | | | | | |
| *** Check Standard: 2 Ck2 | | | | | | | | |
| Line | Flag | %Rcv. | Found | True | Units | SD/RSD | | |
| Ag1 | L | 83.7 | .837 | 1.00 | PPM | .072 | | |
| Ba1 | | 106. | 1.06 | 1.00 | PPM | .010 | | |
| Al3 | | 101. | 1.01 | 1.00 | PPM | .063 | | |
| ----- | | | | | | | | |
| *** Update Slope | | | | | | | | |
| | | | % Corr. | Avg. | SD/RSD | Old | | |
| Ag1 | 1.00 | PPM | 42.9 | .700 | .057 | .700 | | |
| ----- | | | | | | | | |
| *** Check Standard: 2 Ck2 | | | | | | | | |
| Line | Flag | %Rcv. | Found | True | Units | SD/RSD | | |
| Ag1 | | 107. | 1.07 | 1.00 | PPM | .092 | | |
| Ba1 | | 105. | 1.05 | 1.00 | PPM | .008 | | |
| Al3 | | 105. | 1.05 | 1.00 | PPM | .021 | | |

P.E.C. Laboratories
16:26:44 16 Oct 1996

Folder: pripl198
Protocol: priorpl1

Page 22

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|------|-------|-------|--------|---|---|---|---|---|
|------|-------|-------|--------|---|---|---|---|---|

* * Check Standard: 3 Ck3 Seq: 483 16:26:44 16 Oct 1996 ICP

| Line | Flag | %Rcv. | Found | True | Units | SD/RSD |
|------|------|-------|-------|------|-------|--------|
| As2 | | 105. | 1.05 | 1.00 | PPM | .038 |
| Be2 | | 104. | 1.04 | 1.00 | PPM | .007 |
| C1 | H | 112. | 1.12 | 1.00 | PPM | .003 |
| Co1 | H | 112. | 1.12 | 1.00 | PPM | .007 |
| Cr4 | H | 112. | 1.12 | 1.00 | PPM | .011 |
| P1 | H | 113. | 1.13 | 1.00 | PPM | .024 |
| M1 | | 107. | 1.07 | 1.00 | PPM | .006 |
| Sb1 | | 102. | 1.02 | 1.00 | PPM | .058 |
| Tl1 | | 109. | 1.09 | 1.00 | PPM | .071 |
| N3 | H | 112. | 1.12 | 1.00 | PPM | .011 |
| Se1 | | 109. | 1.09 | 1.00 | PPM | .052 |
| V2 | | 108. | 1.08 | 1.00 | PPM | .005 |
| C1 | | 105. | 1.05 | 1.00 | PPM | .006 |
| F2 | H | 110. | 1.10 | 1.00 | PPM | .007 |
| Zn1 | | 109. | 1.09 | 1.00 | PPM | .014 |
| M1 | H | 114. | 1.14 | 1.00 | PPM | .019 |

*** Update Slope Seq: 484 16:30:11 16 Oct 1996 ICP

| | | % Corr. | Avg. | SD/RSD | Old | |
|-----|------|---------|-------|--------|------|------|
| A2 | 3.00 | PPM | -7.52 | 3.24 | .056 | 3.24 |
| Be2 | 3.00 | PPM | -3.06 | 3.09 | .022 | 3.09 |
| Ba1 | 3.00 | PPM | .610 | 2.98 | .019 | 2.98 |
| C1 | 3.00 | PPM | -6.03 | 3.19 | .023 | 3.19 |
| C1 | 3.00 | PPM | -5.59 | 3.18 | .019 | 3.18 |
| Cr4 | 4.50 | PPM | -6.08 | 4.79 | .024 | 4.79 |
| P1 | 3.00 | PPM | -5.76 | 3.18 | .009 | 3.18 |
| M1 | 4.50 | PPM | -4.65 | 4.72 | .027 | 4.72 |
| Sb1 | 3.00 | PPM | -3.71 | 3.12 | .046 | 3.12 |
| Tl1 | 3.00 | PPM | -7.70 | 3.25 | .080 | 3.25 |
| N3 | 4.50 | PPM | -6.89 | 4.83 | .029 | 4.83 |
| Se1 | 3.00 | PPM | -7.69 | 3.25 | .107 | 3.25 |
| V2 | 4.50 | PPM | -4.47 | 4.71 | .033 | 4.71 |
| C1 | 4.50 | PPM | -.319 | 4.51 | .031 | 4.51 |
| F2 | 4.50 | PPM | -5.88 | 4.78 | .027 | 4.78 |
| Zn1 | 4.50 | PPM | -6.61 | 4.82 | .021 | 4.82 |
| A13 | 4.50 | PPM | -2.73 | 4.63 | .061 | 4.63 |
| M2 | 30.0 | PPM | .973 | 29.7 | .225 | 29.7 |
| Ca3 | 30.0 | PPM | -1.99 | 30.6 | .186 | 30.6 |
| Mg4 | 30.0 | PPM | -4.13 | 31.3 | .133 | 31.3 |
| H1 | 30.0 | PPM | -1.44 | 30.4 | .287 | 30.4 |
| Ns1 | 3.00 | PPM | -11.3 | 3.38 | .030 | 3.38 |

** Peak Seq: 485 16:31:58 16 Oct 1996

XAp = 213 YAp = 140 Intensity = 7647395

E.C. Laboratories
16:33:18 16 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 23

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|--------------------------|-------|-------|--------|-------|--------|--------|---|---|
| ***** Update Intercept | | | | | | | | |
| | | | Corr. | Avg. | SD/RSD | Old | | |
| As2 | .000 | PPM | -.004 | .004 | .053 | .009 | | |
| Be2 | .000 | PPM | .021 | -.020 | .001 | .001 | | |
| B1 | .000 | PPM | .027 | -.027 | .000 | -.004 | | |
| C1 | .000 | PPM | .046 | -.043 | .001 | .003 | | |
| Co1 | .000 | PPM | .040 | -.037 | .003 | -.001 | | |
| Cr4 | .000 | PPM | .051 | -.048 | .003 | .001 | | |
| Fe1 | .000 | PPM | .043 | -.041 | .012 | .018 | | |
| Mn1 | .000 | PPM | .037 | -.035 | .001 | .001 | | |
| Sb1 | .000 | PPM | .017 | -.016 | .016 | -.007 | | |
| Tl1 | .000 | PPM | .052 | -.048 | .025 | -.018 | | |
| Ni3 | .000 | PPM | .064 | -.059 | .005 | .008 | | |
| Se1 | .000 | PPM | -.024 | .022 | .049 | -.004 | | |
| V2 | .000 | PPM | .033 | -.032 | .001 | .001 | | |
| Ca1 | .000 | PPM | .046 | -.046 | .001 | .005 | | |
| Fe2 | .000 | PPM | .036 | -.034 | .002 | -.003 | | |
| Zn1 | .000 | PPM | .039 | -.036 | .003 | -.001 | | |
| Al3 | .000 | PPM | -.006 | .006 | .024 | -.005 | | |
| Na2 | .000 | PPM | .073 | -.074 | .192 | -.021 | | |
| Ca3 | .000 | PPM | .234 | -.229 | .002 | -.000 | | |
| Mn4 | .000 | PPM | -.042 | .040 | .028 | -.002 | | |
| K1 | .000 | PPM | -.117 | .115 | .573 | .286 | | |
| Mo1 | .000 | PPM | .006 | -.006 | .006 | -.015 | | |
| >* Check Standard: 2 Ck2 | | | | | | | | |
| Line | Flag | %Rcv. | Found | True | Units | SD/RSD | | |
| Ag1 | H | 115. | 1.15 | 1.00 | PPM | .072 | | |
| Bi1 | | 105. | 1.05 | 1.00 | PPM | .008 | | |
| Al3 | | 97.5 | .975 | 1.00 | PPM | .051 | | |

E.C. Laboratories
16:42:15 16 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 24

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|---------------------------|-------|-------|--------|------|-------|--------|---|---|
| *** Check Standard: 3 Ck3 | | | | | | | | |
| Line | Flag | %Rcv. | Found | True | Units | SD/RSD | | |
| As2 | | 96.9 | .969 | 1.00 | PPM | .047 | | |
| Be2 | | 97.3 | .973 | 1.00 | PPM | .002 | | |
| C1 | | 101. | 1.01 | 1.00 | PPM | .003 | | |
| C1 | | 101. | 1.01 | 1.00 | PPM | .002 | | |
| Cr4 | | 100. | 1.00 | 1.00 | PPM | .002 | | |
| Fe1 | | 98.0 | .980 | 1.00 | PPM | .017 | | |
| Mg1 | | 96.7 | .967 | 1.00 | PPM | .004 | | |
| Sb1 | | 94.8 | .948 | 1.00 | PPM | .044 | | |
| Tl1 | | 104. | 1.04 | 1.00 | PPM | .087 | | |
| Mn3 | | 97.8 | .978 | 1.00 | PPM | .012 | | |
| Se1 | | 94.8 | .948 | 1.00 | PPM | .100 | | |
| V2 | | 97.9 | .979 | 1.00 | PPM | .004 | | |
| Ca1 | | 102. | 1.02 | 1.00 | PPM | .008 | | |
| Po2 | | 98.9 | .989 | 1.00 | PPM | .006 | | |
| Zn1 | | 97.6 | .976 | 1.00 | PPM | .003 | | |
| Mn1 | | 95.3 | .953 | 1.00 | PPM | .018 | | |
| *** Check Standard: 2 Ck2 | | | | | | | | |
| Line | Flag | %Rcv. | Found | True | Units | SD/RSD | | |
| Al1 | | 103. | 1.03 | 1.00 | PPM | .030 | | |
| Ba1 | | 103. | 1.03 | 1.00 | PPM | .009 | | |
| Al3 | | 95.6 | .956 | 1.00 | PPM | .014 | | |
| *** Check Standard: 4 Ck4 | | | | | | | | |
| Line | Flag | %Rcv. | Found | True | Units | SD/RSD | | |
| Na2 | | 99.8 | 49.9 | 50.0 | PPM | .425 | | |
| Ca3 | | 99.1 | 49.6 | 50.0 | PPM | .351 | | |
| Mg4 | | 100. | 10.0 | 10.0 | PPM | .077 | | |
| K1 | | 96.5 | 9.65 | 10.0 | PPM | .255 | | |

E.C. Laboratories
09:23:44 17 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 1

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|-----------|-----------|-------------|--------|---------|-----|----------|-------------|---|
| * | * | Peak | | Seq: | 521 | 09:23:44 | 17 Oct 1996 | |
| XAp = 216 | YAp = 141 | Intensity = | | 7285214 | | | | |
| * | * | Peak | | Seq: | 522 | 09:24:23 | 17 Oct 1996 | |
| XAp = 216 | YAp = 141 | Intensity = | | 7321908 | | | | |

E.C. Laboratories
09:40:18 17 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 2

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|------------------------|-------|-------|-------------|--------|----------------|--------|--------|---|
| * * Standard: 1 Rep: 1 | | | | | | | | |
| As2 | .000 | PPM | | 965 | 15 | 653 | -723 | |
| | | | Ave. Int. = | | 228 S. D. = | 747 | | |
| B2 | .000 | PPM | | -1585 | -4513 | -3393 | -953 | |
| | | | Ave. Int. = | | -2611 S. D. = | 1636 | | |
| Ba1 | .000 | PPM | | 103254 | 104309 | 101614 | 104972 | |
| | | | Ave. Int. = | | 103537 S. D. = | 1464 | | |
| Cd1 | .000 | PPM | | -5235 | -2777 | 893 | -1875 | |
| | | | Ave. Int. = | | -2249 S. D. = | 2530 | | |
| Cr1 | .000 | PPM | | -505 | 583 | -433 | 927 | |
| | | | Ave. Int. = | | 143 S. D. = | 721 | | |
| Cr4 | .000 | PPM | | -4063 | -5951 | -6825 | -2903 | |
| | | | Ave. Int. = | | -4936 S. D. = | 1779 | | |
| F1 | .000 | PPM | | -899 | -1643 | -2185 | 1847 | |
| | | | Ave. Int. = | | -720 S. D. = | 1791 | | |
| Mn1 | .000 | PPM | | 13593 | 12961 | 11289 | 8497 | |
| | | | Ave. Int. = | | 11585 S. D. = | 2277 | | |
| S1 | .000 | PPM | | 767 | 129 | -297 | -295 | |
| | | | Ave. Int. = | | 76 S. D. = | 502 | | |
| T1 | .000 | PPM | | -533 | -101 | 223 | 63 | |
| | | | Ave. Int. = | | -87 S. D. = | 325 | | |
| Ni3 | .000 | PPM | | 2805 | 565 | 2359 | 2815 | |
| | | | Ave. Int. = | | 2136 S. D. = | 1069 | | |
| S1 | .000 | PPM | | 703 | 413 | 583 | 823 | |
| | | | Ave. Int. = | | 631 S. D. = | 175 | | |
| V2 | .000 | PPM | | 147929 | 154739 | 157825 | 154281 | |
| | | | Ave. Int. = | | 153694 S. D. = | 4153 | | |
| Cu1 | .000 | PPM | | 267 | 5257 | 5897 | 2785 | |
| | | | Ave. Int. = | | 3552 S. D. = | 2568 | | |
| Fe2 | .000 | PPM | | 26791 | 22711 | 30577 | 23247 | |
| | | | Ave. Int. = | | 25832 S. D. = | 3645 | | |
| Zn1 | .000 | PPM | | 4637 | 4277 | 2725 | 4439 | |
| | | | Ave. Int. = | | 4020 S. D. = | 875 | | |
| Al3 | .000 | PPM | | 11911 | 10567 | 13345 | 13743 | |
| | | | Ave. Int. = | | 12392 S. D. = | 1449 | | |
| Na2 | .000 | PPM | | -4611 | -2691 | -4051 | -3435 | |
| | | | Ave. Int. = | | -3697 S. D. = | 825 | | |
| Ca3 | .000 | PPM | | 33830 | 41428 | 45196 | 38662 | |
| | | | Ave. Int. = | | 39779 S. D. = | 4785 | | |
| Mg4 | .000 | PPM | | 22132 | 21150 | 23916 | 21790 | |
| | | | Ave. Int. = | | 22247 S. D. = | 1185 | | |
| K1 | .000 | PPM | | -1083 | -1183 | -1505 | -977 | |
| | | | Ave. Int. = | | -1187 S. D. = | 228 | | |
| Li1 | .000 | PPM | | 1279 | 1573 | 1071 | 1567 | |
| | | | Ave. Int. = | | 1373 S. D. = | 243 | | |

E.C. Laboratories
09:41:13 17 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 3

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|----------------------|-------|-------|-------------|--------|----------------|--------|--------|---|
| * Standard: 1 Rep: 2 | | | | | | | | |
| As2 | .000 | PPM | | 1181 | -436 | 827 | -1227 | |
| | | | Ave. Int. = | | 86 S. D. = | | 1117 | |
| B2 | .000 | PPM | | -3081 | -3569 | -891 | -179 | |
| | | | Ave. Int. = | | -1930 S. D. = | | 1649 | |
| Ba1 | .000 | PPM | | 115312 | 113598 | 102148 | 106164 | |
| | | | Ave. Int. = | | 109306 S. D. = | | 6208 | |
| Cd1 | .000 | PPM | | 335 | -583 | -603 | 1223 | |
| | | | Ave. Int. = | | 93 S. D. = | | 871 | |
| Cr1 | .000 | PPM | | 711 | -145 | 837 | -1292 | |
| | | | Ave. Int. = | | 28 S. D. = | | 982 | |
| Cr4 | .000 | PPM | | -703 | -9343 | -6527 | -1015 | |
| | | | Ave. Int. = | | -4397 S. D. = | | 4246 | |
| Fe1 | .000 | PPM | | -705 | -3377 | -729 | -353 | |
| | | | Ave. Int. = | | -1291 S. D. = | | 1401 | |
| Mn1 | .000 | PPM | | 9601 | 7857 | 9921 | 15105 | |
| | | | Ave. Int. = | | 10621 S. D. = | | 3124 | |
| Si1 | .000 | PPM | | 1211 | 675 | 471 | 99 | |
| | | | Ave. Int. = | | 614 S. D. = | | 464 | |
| Tl1 | .000 | PPM | | -9 | -377 | 35 | -141 | |
| | | | Ave. Int. = | | -123 S. D. = | | 185 | |
| Ni3 | .000 | PPM | | 2703 | 2749 | -585 | 3535 | |
| | | | Ave. Int. = | | 2101 S. D. = | | 1831 | |
| Se1 | .000 | PPM | | 317 | -209 | 751 | 375 | |
| | | | Ave. Int. = | | 309 S. D. = | | 395 | |
| V2 | .000 | PPM | | 157241 | 164203 | 159185 | 160745 | |
| | | | Ave. Int. = | | 160344 S. D. = | | 2945 | |
| Al1 | .000 | PPM | | 5285 | 6235 | -3689 | 2999 | |
| | | | Ave. Int. = | | 2708 S. D. = | | 4475 | |
| Fe2 | .000 | PPM | | 21217 | 20033 | 18753 | 35665 | |
| | | | Ave. Int. = | | 23917 S. D. = | | 7896 | |
| Zn1 | .000 | PPM | | 4173 | 3869 | 4165 | 4845 | |
| | | | Ave. Int. = | | 4263 S. D. = | | 413 | |
| Li3 | .000 | PPM | | 12229 | 4657 | 12479 | 5649 | |
| | | | Ave. Int. = | | 8754 S. D. = | | 4178 | |
| Na2 | .000 | PPM | | -4521 | -3267 | -1955 | -3791 | |
| | | | Ave. Int. = | | -3384 S. D. = | | 1082 | |
| Ca3 | .000 | PPM | | 58630 | 43886 | 40372 | 36958 | |
| | | | Ave. Int. = | | 44962 S. D. = | | 9541 | |
| Mg4 | .000 | PPM | | 23836 | 21836 | 22044 | 20966 | |
| | | | Ave. Int. = | | 22171 S. D. = | | 1205 | |
| K1 | .000 | PPM | | -881 | -1057 | -1101 | -1361 | |
| | | | Ave. Int. = | | -1100 S. D. = | | 198 | |
| Ca1 | .000 | PPM | | 1447 | 2047 | 2373 | 2437 | |
| | | | Ave. Int. = | | 2076 S. D. = | | 453 | |

E.C. Laboratories
09:42:09 17 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 4

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|----------------------|-------|-------|-------------|--------|----------------|--------|--------|---|
| * Standard: 1 Rep: 3 | | | | | | | | |
| As2 | .000 | PPM | | -851 | -969 | 543 | 501 | |
| | | | Ave. Int. = | | -194 S. D. = | 828 | | |
| B2 | .000 | PPM | | -5163 | -2859 | -2113 | -2947 | |
| | | | Ave. Int. = | | -3271 S. D. = | 1316 | | |
| Ba1 | .000 | PPM | | 99614 | 104564 | 108248 | 110046 | |
| | | | Ave. Int. = | | 105618 S. D. = | 4607 | | |
| Cd1 | .000 | PPM | | 263 | 487 | -305 | 2551 | |
| | | | Ave. Int. = | | 749 S. D. = | 1247 | | |
| Co1 | .000 | PPM | | 637 | -9 | 159 | -257 | |
| | | | Ave. Int. = | | 133 S. D. = | 377 | | |
| Cr4 | .000 | PPM | | 951 | -6473 | -12855 | -1337 | |
| | | | Ave. Int. = | | -4929 S. D. = | 6129 | | |
| Fe1 | .000 | PPM | | -755 | -2209 | 405 | -2521 | |
| | | | Ave. Int. = | | -1270 S. D. = | 1356 | | |
| Mn1 | .000 | PPM | | 7737 | 5833 | 6827 | 6434 | |
| | | | Ave. Int. = | | 6708 S. D. = | 799 | | |
| Se1 | .000 | PPM | | 959 | 263 | -245 | -323 | |
| | | | Ave. Int. = | | 164 S. D. = | 591 | | |
| Tl1 | .000 | PPM | | -639 | 3 | -505 | -319 | |
| | | | Ave. Int. = | | -365 S. D. = | 278 | | |
| Ni3 | .000 | PPM | | 3549 | 2535 | 1191 | 3629 | |
| | | | Ave. Int. = | | 2726 S. D. = | 1138 | | |
| Si1 | .000 | PPM | | 495 | 1173 | -1 | 487 | |
| | | | Ave. Int. = | | 539 S. D. = | 482 | | |
| V2 | .000 | PPM | | 157035 | 157075 | 160169 | 161913 | |
| | | | Ave. Int. = | | 159048 S. D. = | 2409 | | |
| Cu1 | .000 | PPM | | 6211 | 2983 | 4081 | 4713 | |
| | | | Ave. Int. = | | 4497 S. D. = | 1348 | | |
| Fe2 | .000 | PPM | | 26753 | 19113 | 18793 | 18199 | |
| | | | Ave. Int. = | | 20715 S. D. = | 4043 | | |
| Zn1 | .000 | PPM | | 3445 | 4661 | 2911 | 2621 | |
| | | | Ave. Int. = | | 3410 S. D. = | 901 | | |
| Al3 | .000 | PPM | | 13697 | 10007 | 11871 | 11943 | |
| | | | Ave. Int. = | | 11880 S. D. = | 1507 | | |
| Na2 | .000 | PPM | | -4691 | -4401 | -1867 | -3787 | |
| | | | Ave. Int. = | | -3687 S. D. = | 1270 | | |
| Ca3 | .000 | PPM | | 36782 | 42094 | 49468 | 31522 | |
| | | | Ave. Int. = | | 39967 S. D. = | 7665 | | |
| Mg4 | .000 | PPM | | 22646 | 20078 | 22446 | 21374 | |
| | | | Ave. Int. = | | 21636 S. D. = | 1179 | | |
| K1 | .000 | PPM | | -305 | -587 | -745 | -1477 | |
| | | | Ave. Int. = | | -779 S. D. = | 500 | | |
| Li1 | .000 | PPM | | 1039 | 39 | 2199 | 1965 | |
| | | | Ave. Int. = | | 1311 S. D. = | 985 | | |

E.C. Laboratories
09:45:18 17 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 5

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|------------------------|-------|-------|-------------|----------|----------|----------|----------|---|
| * * Standard: 2 Rep: 1 | | | | | | | | |
| As2 | 1.00 | PPM | | 13077 | 15399 | 13821 | 14783 | |
| | | | Ave. Int. = | 14270 | S. D. = | 1027 | | |
| B2 | 1.00 | PPM | | 3801317 | 3916468 | 3871898 | 3863896 | |
| | | | Ave. Int. = | 3863395 | S. D. = | 47409 | | |
| Ba1 | 1.00 | PPM | | 8205612 | 8464003 | 8371204 | 8383817 | |
| | | | Ave. Int. = | 8356159 | S. D. = | 108453 | | |
| Cu1 | 1.00 | PPM | | 1644442 | 1688422 | 1677964 | 1662598 | |
| | | | Ave. Int. = | 1668357 | S. D. = | 19148 | | |
| Co1 | 1.00 | PPM | | 184012 | 189078 | 186538 | 187498 | |
| | | | Ave. Int. = | 186782 | S. D. = | 2123 | | |
| Cr4 | 1.50 | PPM | | 2180939 | 2231328 | 2207881 | 2203403 | |
| | | | Ave. Int. = | 2205888 | S. D. = | 20654 | | |
| Fe1 | 1.00 | PPM | | 57189 | 55557 | 57943 | 55077 | |
| | | | Ave. Int. = | 56442 | S. D. = | 1349 | | |
| Mn1 | 1.50 | PPM | | 3654284 | 3753727 | 3714290 | 3708196 | |
| | | | Ave. Int. = | 3707624 | S. D. = | 40888 | | |
| Si1 | 1.00 | PPM | | 18559 | 19205 | 19221 | 18433 | |
| | | | Ave. Int. = | 18855 | S. D. = | 417 | | |
| Tl1 | 1.00 | PPM | | 4013 | 4831 | 4535 | 4399 | |
| | | | Ave. Int. = | 4445 | S. D. = | 340 | | |
| Ni3 | 1.50 | PPM | | 357621 | 367669 | 369011 | 362741 | |
| | | | Ave. Int. = | 364261 | S. D. = | 5183 | | |
| Se1 | 1.00 | PPM | | 7583 | 7711 | 7493 | 6885 | |
| | | | Ave. Int. = | 7418 | S. D. = | 366 | | |
| V2 | 1.50 | PPM | | 2824006 | 2886994 | 2871434 | 2866970 | |
| | | | Ave. Int. = | 2862351 | S. D. = | 26966 | | |
| Ca1 | 1.50 | PPM | | 1081538 | 1106160 | 1099302 | 1103664 | |
| | | | Ave. Int. = | 1097666 | S. D. = | 11119 | | |
| Fe2 | 1.50 | PPM | | 1121826 | 1155612 | 1146386 | 1145436 | |
| | | | Ave. Int. = | 1142315 | S. D. = | 14410 | | |
| Zn1 | 1.50 | PPM | | 490120 | 496704 | 496320 | 496778 | |
| | | | Ave. Int. = | 494981 | S. D. = | 3247 | | |
| Al3 | 1.50 | PPM | | 149285 | 152831 | 151159 | 148455 | |
| | | | Ave. Int. = | 150433 | S. D. = | 1959 | | |
| Na2 | 10.0 | PPM | | 111392 | 115304 | 116120 | 115734 | |
| | | | Ave. Int. = | 114638 | S. D. = | 2189 | | |
| Ca3 | 10.0 | PPM | | 10326444 | 10605970 | 10503219 | 10505739 | |
| | | | Ave. Int. = | 10485343 | S. D. = | 116240 | | |
| Mg4 | 10.0 | PPM | | 495076 | 511268 | 500634 | 504378 | |
| | | | Ave. Int. = | 502839 | S. D. = | 6796 | | |
| K1 | 10.0 | PPM | | 9148 | 8450 | 8622 | 8584 | |
| | | | Ave. Int. = | 8701 | S. D. = | 307 | | |
| Cl1 | 1.00 | PPM | | 79272 | 80264 | 78272 | 77546 | |
| | | | Ave. Int. = | 78839 | S. D. = | 1185 | | |

E.C. Laboratories
09:46:15 17 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 6

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|------|---------------|--------|-------------|----------|--------------------------|----------|----------|---|
| * | * Standard: 2 | Rep: 2 | | Seq: 527 | 09:46:15 17 Oct 1996 ICP | | | |
| As2 | 1.00 | PPM | | 14469 | 15669 | 13783 | 14341 | |
| | | | Ave. Int. = | 14566 | S. D. = | 794 | | |
| B2 | 1.00 | PPM | | 3850648 | 3831950 | 3766890 | 3795307 | |
| | | | Ave. Int. = | 3811199 | S. D. = | 37428 | | |
| Ba1 | 1.00 | PPM | | 8342388 | 8296449 | 8118342 | 8187481 | |
| | | | Ave. Int. = | 8236165 | S. D. = | 101931 | | |
| Cd1 | 1.00 | PPM | | 1662582 | 1660812 | 1635752 | 1645032 | |
| | | | Ave. Int. = | 1651045 | S. D. = | 12891 | | |
| Cl1 | 1.00 | PPM | | 187282 | 187980 | 183540 | 183484 | |
| | | | Ave. Int. = | 185572 | S. D. = | 2395 | | |
| Cr4 | 1.50 | PPM | | 2190421 | 2207625 | 2148699 | 2159611 | |
| | | | Ave. Int. = | 2176589 | S. D. = | 27207 | | |
| P1 | 1.00 | PPM | | 53885 | 56847 | 56695 | 55669 | |
| | | | Ave. Int. = | 55774 | S. D. = | 1364 | | |
| Mn1 | 1.50 | PPM | | 3691394 | 3684644 | 3618925 | 3636317 | |
| | | | Ave. Int. = | 3657820 | S. D. = | 35693 | | |
| S1 | 1.00 | PPM | | 19363 | 18711 | 18453 | 16471 | |
| | | | Ave. Int. = | 18250 | S. D. = | 1246 | | |
| T1 | 1.00 | PPM | | 4119 | 4211 | 4573 | 4693 | |
| | | | Ave. Int. = | 4399 | S. D. = | 277 | | |
| Ni3 | 1.50 | PPM | | 362603 | 362389 | 357013 | 361885 | |
| | | | Ave. Int. = | 360973 | S. D. = | 2657 | | |
| S1 | 1.00 | PPM | | 6831 | 7741 | 7357 | 7125 | |
| | | | Ave. Int. = | 7264 | S. D. = | 384 | | |
| V2 | 1.50 | PPM | | 2855066 | 2864946 | 2795950 | 2816144 | |
| | | | Ave. Int. = | 2833027 | S. D. = | 32477 | | |
| Ca1 | 1.50 | PPM | | 1093852 | 1087386 | 1071412 | 1082942 | |
| | | | Ave. Int. = | 1083898 | S. D. = | 9453 | | |
| F2 | 1.50 | PPM | | 1139954 | 1140012 | 1115652 | 1124820 | |
| | | | Ave. Int. = | 1130110 | S. D. = | 12000 | | |
| Zn1 | 1.50 | PPM | | 496122 | 495474 | 482290 | 483986 | |
| | | | Ave. Int. = | 489468 | S. D. = | 7347 | | |
| Al3 | 1.50 | PPM | | 150869 | 149559 | 147725 | 145437 | |
| | | | Ave. Int. = | 148398 | S. D. = | 2358 | | |
| Na2 | 10.0 | PPM | | 114760 | 114760 | 115654 | 113464 | |
| | | | Ave. Int. = | 114660 | S. D. = | 902 | | |
| Ca3 | 10.0 | PPM | | 10458386 | 10454843 | 10222232 | 10295290 | |
| | | | Ave. Int. = | 10357688 | S. D. = | 118069 | | |
| Mg4 | 10.0 | PPM | | 497660 | 501194 | 490364 | 489260 | |
| | | | Ave. Int. = | 494620 | S. D. = | 5753 | | |
| K1 | 10.0 | PPM | | 7866 | 8370 | 7926 | 7430 | |
| | | | Ave. Int. = | 7898 | S. D. = | 385 | | |
| Na1 | 1.00 | PPM | | 79504 | 76968 | 76792 | 78046 | |
| | | | Ave. Int. = | 77828 | S. D. = | 1248 | | |

E.C. Laboratories
09:47:12 17 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 7

| Name | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|------------------------|-------|-------|-------------|----------|--------------------------|----------|----------|---|
| * * Standard: 2 Rep: 3 | | | | Seq: 528 | 09:47:12 17 Oct 1996 ICP | | | |
| As2 | 1.00 | PPM | | 14055 | 15535 | 14101 | 14341 | |
| | | | Ave. Int. = | 14508 | S. D. = | 696 | | |
| Se2 | 1.00 | PPM | | 3823459 | 3817035 | 3726946 | 3788451 | |
| | | | Ave. Int. = | 3788973 | S. D. = | 44062 | | |
| Ba1 | 1.00 | PPM | | 8269567 | 8269785 | 8039308 | 8197427 | |
| | | | Ave. Int. = | 8194022 | S. D. = | 108620 | | |
| Ca1 | 1.00 | PPM | | 1645896 | 1648546 | 1622058 | 1646792 | |
| | | | Ave. Int. = | 1640823 | S. D. = | 12558 | | |
| Cl1 | 1.00 | PPM | | 186450 | 185236 | 180178 | 185524 | |
| | | | Ave. Int. = | 184347 | S. D. = | 2827 | | |
| Cr4 | 1.50 | PPM | | 2186675 | 2184985 | 2140307 | 2171569 | |
| | | | Ave. Int. = | 2170884 | S. D. = | 21476 | | |
| Fe1 | 1.00 | PPM | | 55941 | 54911 | 51799 | 57479 | |
| | | | Ave. Int. = | 55033 | S. D. = | 2400 | | |
| Mn1 | 1.50 | PPM | | 3683236 | 3671468 | 3584037 | 3636141 | |
| | | | Ave. Int. = | 3643721 | S. D. = | 44538 | | |
| Si1 | 1.00 | PPM | | 18493 | 17253 | 17899 | 19421 | |
| | | | Ave. Int. = | 18267 | S. D. = | 921 | | |
| Tl1 | 1.00 | PPM | | 3639 | 4415 | 4169 | 4095 | |
| | | | Ave. Int. = | 4080 | S. D. = | 324 | | |
| Ni3 | 1.50 | PPM | | 359957 | 360139 | 351133 | 359387 | |
| | | | Ave. Int. = | 357654 | S. D. = | 4359 | | |
| Sn1 | 1.00 | PPM | | 7199 | 7301 | 8119 | 6383 | |
| | | | Ave. Int. = | 7251 | S. D. = | 710 | | |
| V2 | 1.50 | PPM | | 2829914 | 2830090 | 2776118 | 2818824 | |
| | | | Ave. Int. = | 2813737 | S. D. = | 25627 | | |
| Ca1 | 1.50 | PPM | | 1088752 | 1088300 | 1059530 | 1078604 | |
| | | | Ave. Int. = | 1078797 | S. D. = | 13671 | | |
| Fe2 | 1.50 | PPM | | 1142914 | 1134330 | 1103834 | 1120730 | |
| | | | Ave. Int. = | 1125452 | S. D. = | 17062 | | |
| Zn1 | 1.50 | PPM | | 490240 | 487216 | 478948 | 491634 | |
| | | | Ave. Int. = | 487010 | S. D. = | 5682 | | |
| Al3 | 1.50 | PPM | | 148749 | 145279 | 145269 | 151871 | |
| | | | Ave. Int. = | 147792 | S. D. = | 3175 | | |
| Na2 | 10.0 | PPM | | 114224 | 114008 | 111080 | 111960 | |
| | | | Ave. Int. = | 112818 | S. D. = | 1544 | | |
| Ca3 | 10.0 | PPM | | 10410314 | 10376653 | 10154881 | 10303514 | |
| | | | Ave. Int. = | 10311341 | S. D. = | 113435 | | |
| Mg4 | 10.0 | PPM | | 497500 | 497908 | 488324 | 498828 | |
| | | | Ave. Int. = | 495640 | S. D. = | 4909 | | |
| K1 | 10.0 | PPM | | 8386 | 8806 | 8562 | 8322 | |
| | | | Ave. Int. = | 8519 | S. D. = | 217 | | |
| Mo1 | 1.00 | PPM | | 79350 | 77512 | 75478 | 77336 | |
| | | | Ave. Int. = | 77419 | S. D. = | 1582 | | |

E.C. Laboratories
09:49:01 17 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 8

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|------|-------|-------|--------|---|---|---|---|---|
|------|-------|-------|--------|---|---|---|---|---|

* Peak Seq: 529 09:49:01 17 Oct 1996

XAp = 215 YAp = 141 Intensity = 7464047

S.E.C. Laboratories
09:50:21 17 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 9

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|-----------------------|-------|-------|-------------|----------|----------|----------|----------|---|
| >* Standard: 3 Rep: 1 | | | | | | | | |
| As2 | 3.00 | PPM | | 44367 | 43159 | 42597 | 45431 | |
| | | | Ave. Int. = | 43889 | S. D. = | 1266 | | |
| B2 | 3.00 | PPM | | 11379951 | 11315789 | 11396630 | 11603171 | |
| | | | Ave. Int. = | 11423885 | S. D. = | 124501 | | |
| Ba1 | 3.00 | PPM | | 24665755 | 24501866 | 24693163 | 25153632 | |
| | | | Ave. Int. = | 24753604 | S. D. = | 279741 | | |
| Cd1 | 3.00 | PPM | | 4934601 | 4913352 | 4948575 | 5023602 | |
| | | | Ave. Int. = | 4955033 | S. D. = | 47952 | | |
| Co1 | 3.00 | PPM | | 549385 | 544185 | 549103 | 558129 | |
| | | | Ave. Int. = | 550201 | S. D. = | 5800 | | |
| Cr4 | 4.50 | PPM | | 6429995 | 6394776 | 6423749 | 6549934 | |
| | | | Ave. Int. = | 6449614 | S. D. = | 68618 | | |
| Fe1 | 3.00 | PPM | | 166213 | 162923 | 165037 | 167861 | |
| | | | Ave. Int. = | 165509 | S. D. = | 2077 | | |
| Mn1 | 4.50 | PPM | | 10870628 | 10821925 | 10871608 | 11053239 | |
| | | | Ave. Int. = | 10904350 | S. D. = | 101933 | | |
| Mo1 | 3.00 | PPM | | 55859 | 55637 | 54751 | 56663 | |
| | | | Ave. Int. = | 55728 | S. D. = | 786 | | |
| Tl1 | 3.00 | PPM | | 13103 | 13275 | 13053 | 13891 | |
| | | | Ave. Int. = | 13331 | S. D. = | 386 | | |
| Ni3 | 4.50 | PPM | | 1048880 | 1037138 | 1048912 | 1068960 | |
| | | | Ave. Int. = | 1050973 | S. D. = | 13211 | | |
| Se1 | 3.00 | PPM | | 20709 | 19901 | 20511 | 20135 | |
| | | | Ave. Int. = | 20314 | S. D. = | 364 | | |
| V2 | 4.50 | PPM | | 7998802 | 7964986 | 8009602 | 8171501 | |
| | | | Ave. Int. = | 8036223 | S. D. = | 92166 | | |
| Al1 | 4.50 | PPM | | 3242217 | 3235625 | 3254525 | 3310700 | |
| | | | Ave. Int. = | 3260767 | S. D. = | 34198 | | |
| Fe2 | 4.50 | PPM | | 3284867 | 3279435 | 3289905 | 3352664 | |
| | | | Ave. Int. = | 3301718 | S. D. = | 34232 | | |
| Zn1 | 4.50 | PPM | | 1432679 | 1421821 | 1435755 | 1461149 | |
| | | | Ave. Int. = | 1437851 | S. D. = | 16642 | | |
| Li3 | 4.50 | PPM | | 412695 | 422279 | 422071 | 431119 | |
| | | | Ave. Int. = | 422041 | S. D. = | 7524 | | |
| Na2 | 30.0 | PPM | | 348176 | 347950 | 347080 | 358304 | |
| | | | Ave. Int. = | 350378 | S. D. = | 5305 | | |
| La3 | 30.0 | PPM | | 31102078 | 30974819 | 31209191 | 31729648 | |
| | | | Ave. Int. = | 31253934 | S. D. = | 331296 | | |
| Mg4 | 30.0 | PPM | | 1433524 | 1431348 | 1440282 | 1450356 | |
| | | | Ave. Int. = | 1438878 | S. D. = | 8546 | | |
| K1 | 30.0 | PPM | | 26777 | 27915 | 27305 | 27093 | |
| | | | Ave. Int. = | 27273 | S. D. = | 480 | | |
| Co1 | 3.00 | PPM | | 217065 | 217011 | 217401 | 225209 | |
| | | | Ave. Int. = | 219172 | S. D. = | 4029 | | |

E.C. Laboratories
09:51:15 17 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 10

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|------------------------|-------|-------|-------------|----------|--------------------------|----------|----------|---|
| * * Standard: 3 Rep: 2 | | | | | | | | |
| | | | | Seq: 531 | 09:51:15 17 Oct 1996 ICP | | | |
| As2 | 3.00 | PPM | | 43421 | 44445 | 43741 | 44393 | |
| | | | Ave. Int. = | 44000 | S. D. = | 502 | | |
| Ba2 | 3.00 | PPM | | 11397600 | 11352272 | 11404780 | 11408310 | |
| | | | Ave. Int. = | 11390741 | S. D. = | 26030 | | |
| Ba1 | 3.00 | PPM | | 24696851 | 24570895 | 24701571 | 24703197 | |
| | | | Ave. Int. = | 24668129 | S. D. = | 64878 | | |
| Cd1 | 3.00 | PPM | | 4958759 | 4936169 | 4976881 | 4963503 | |
| | | | Ave. Int. = | 4958828 | S. D. = | 16943 | | |
| Co1 | 3.00 | PPM | | 552921 | 545449 | 551975 | 551201 | |
| | | | Ave. Int. = | 550387 | S. D. = | 3366 | | |
| Cr4 | 4.50 | PPM | | 6431147 | 6406829 | 6434427 | 6452547 | |
| | | | Ave. Int. = | 6431238 | S. D. = | 18798 | | |
| Fe1 | 3.00 | PPM | | 164997 | 159885 | 167029 | 170811 | |
| | | | Ave. Int. = | 165681 | S. D. = | 4553 | | |
| Mn1 | 4.50 | PPM | | 10855748 | 10853428 | 10889476 | 10888694 | |
| | | | Ave. Int. = | 10871837 | S. D. = | 19942 | | |
| Se1 | 3.00 | PPM | | 56099 | 53365 | 54769 | 55679 | |
| | | | Ave. Int. = | 54978 | S. D. = | 1210 | | |
| Tl1 | 3.00 | PPM | | 13995 | 13103 | 13511 | 13397 | |
| | | | Ave. Int. = | 13502 | S. D. = | 371 | | |
| Ni3 | 4.50 | PPM | | 1052552 | 1045994 | 1049904 | 1053552 | |
| | | | Ave. Int. = | 1050501 | S. D. = | 3376 | | |
| Se1 | 3.00 | PPM | | 19797 | 20583 | 20207 | 20821 | |
| | | | Ave. Int. = | 20352 | S. D. = | 448 | | |
| V2 | 4.50 | PPM | | 8016092 | 7981794 | 8027785 | 8036333 | |
| | | | Ave. Int. = | 8015501 | S. D. = | 23954 | | |
| U1 | 4.50 | PPM | | 3259661 | 3235341 | 3255465 | 3249373 | |
| | | | Ave. Int. = | 3249960 | S. D. = | 10622 | | |
| Fe2 | 4.50 | PPM | | 3291211 | 3276067 | 3282379 | 3301169 | |
| | | | Ave. Int. = | 3287707 | S. D. = | 10915 | | |
| Zn1 | 4.50 | PPM | | 1440853 | 1433437 | 1441751 | 1437303 | |
| | | | Ave. Int. = | 1438336 | S. D. = | 3789 | | |
| Li3 | 4.50 | PPM | | 421157 | 415725 | 419229 | 427037 | |
| | | | Ave. Int. = | 420787 | S. D. = | 4735 | | |
| Na2 | 30.0 | PPM | | 352360 | 349936 | 349921 | 347416 | |
| | | | Ave. Int. = | 349908 | S. D. = | 2019 | | |
| La3 | 30.0 | PPM | | 31171911 | 31014086 | 31228599 | 31202265 | |
| | | | Ave. Int. = | 31154215 | S. D. = | 96248 | | |
| Mg4 | 30.0 | PPM | | 1444330 | 1434266 | 1439548 | 1437666 | |
| | | | Ave. Int. = | 1438953 | S. D. = | 4199 | | |
| K1 | 30.0 | PPM | | 27653 | 26873 | 28123 | 27133 | |
| | | | Ave. Int. = | 27446 | S. D. = | 556 | | |
| Ca1 | 3.00 | PPM | | 219469 | 217081 | 220361 | 220849 | |
| | | | Ave. Int. = | 219440 | S. D. = | 1673 | | |

E.C. Laboratories
09:52:12 17 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 11

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|------|---------------|--------|-------------|----------|--------------------------|----------|----------|---|
| * | * Standard: 3 | Rep: 3 | | Seq: 532 | 09:52:12 17 Oct 1996 ICP | | | |
| As2 | 3.00 | PPM | | 46511 | 44669 | 44921 | 44783 | |
| | | | Ave. Int. = | 45221 | S. D. = | 866 | | |
| B2 | 3.00 | PPM | | 11756775 | 11457278 | 11738920 | 11527892 | |
| | | | Ave. Int. = | 11620216 | S. D. = | 150346 | | |
| Ba1 | 3.00 | PPM | | 25555443 | 24884839 | 25525416 | 24967424 | |
| | | | Ave. Int. = | 25233281 | S. D. = | 356475 | | |
| Ca1 | 3.00 | PPM | | 5059337 | 4958065 | 5048842 | 4970081 | |
| | | | Ave. Int. = | 5009081 | S. D. = | 52378 | | |
| Cl1 | 3.00 | PPM | | 564751 | 553425 | 564081 | 552449 | |
| | | | Ave. Int. = | 558677 | S. D. = | 6645 | | |
| Cr4 | 4.50 | PPM | | 6626670 | 6458773 | 6607174 | 6479299 | |
| | | | Ave. Int. = | 6542979 | S. D. = | 86161 | | |
| Fe1 | 3.00 | PPM | | 167637 | 162235 | 165811 | 165317 | |
| | | | Ave. Int. = | 165250 | S. D. = | 2244 | | |
| Mn1 | 4.50 | PPM | | 11198937 | 10920990 | 11175680 | 10972972 | |
| | | | Ave. Int. = | 11067145 | S. D. = | 140687 | | |
| S1 | 3.00 | PPM | | 59333 | 55365 | 55233 | 55089 | |
| | | | Ave. Int. = | 56255 | S. D. = | 2055 | | |
| Tl1 | 3.00 | PPM | | 14280 | 14401 | 13647 | 13623 | |
| | | | Ave. Int. = | 13988 | S. D. = | 410 | | |
| Ni3 | 4.50 | PPM | | 1075200 | 1052944 | 1070280 | 1053186 | |
| | | | Ave. Int. = | 1062903 | S. D. = | 11536 | | |
| S1 | 3.00 | PPM | | 21359 | 19655 | 21541 | 19461 | |
| | | | Ave. Int. = | 20504 | S. D. = | 1098 | | |
| V2 | 4.50 | PPM | | 8252097 | 8078829 | 8241055 | 8099192 | |
| | | | Ave. Int. = | 8167793 | S. D. = | 91461 | | |
| Ca1 | 4.50 | PPM | | 3358118 | 3283759 | 3350162 | 3291847 | |
| | | | Ave. Int. = | 3320972 | S. D. = | 38579 | | |
| Fe2 | 4.50 | PPM | | 3395696 | 3314713 | 3398384 | 3324154 | |
| | | | Ave. Int. = | 3358237 | S. D. = | 44985 | | |
| Zn1 | 4.50 | PPM | | 1479880 | 1443335 | 1471021 | 1444935 | |
| | | | Ave. Int. = | 1459793 | S. D. = | 18450 | | |
| Al3 | 4.50 | PPM | | 431725 | 424191 | 433727 | 426031 | |
| | | | Ave. Int. = | 428919 | S. D. = | 4534 | | |
| Na2 | 30.0 | PPM | | 365128 | 352998 | 359684 | 354190 | |
| | | | Ave. Int. = | 358000 | S. D. = | 5573 | | |
| Ca3 | 30.0 | PPM | | 32057613 | 31283547 | 31964971 | 31337285 | |
| | | | Ave. Int. = | 31660854 | S. D. = | 407006 | | |
| Mg4 | 30.0 | PPM | | 1479591 | 1446002 | 1470871 | 1447066 | |
| | | | Ave. Int. = | 1460883 | S. D. = | 16952 | | |
| K1 | 30.0 | PPM | | 28265 | 27407 | 28795 | 28229 | |
| | | | Ave. Int. = | 28174 | S. D. = | 573 | | |
| Li1 | 3.00 | PPM | | 228817 | 221449 | 225291 | 216714 | |
| | | | Ave. Int. = | 223068 | S. D. = | 5196 | | |

E.C. Laboratories
09:55:21 17 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 12

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|----------------------|-------|-------|-------------|----------|--------------------------|----------|----------|---|
| * Standard: 4 Rep: 1 | | | | | | | | |
| | | | | Seq: 533 | 09:55:21 17 Oct 1996 ICP | | | |
| As2 | 7.00 | PPM | | 104410 | 103900 | 103142 | 104530 | |
| | | | Ave. Int. = | 103996 | S. D. = | 631 | | |
| Ba2 | 7.00 | PPM | | 26635704 | 26855206 | 26626120 | 26584534 | |
| | | | Ave. Int. = | 26675391 | S. D. = | 121917 | | |
| Ba1 | 7.00 | PPM | | 57577926 | 58040358 | 57413597 | 57345002 | |
| | | | Ave. Int. = | 57594221 | S. D. = | 313070 | | |
| Cd1 | 7.00 | PPM | | 11324778 | 11427980 | 11288985 | 11284419 | |
| | | | Ave. Int. = | 11331541 | S. D. = | 66778 | | |
| Cr1 | 7.00 | PPM | | 1266267 | 1281585 | 1267397 | 1263049 | |
| | | | Ave. Int. = | 1269575 | S. D. = | 8216 | | |
| Cr4 | 10.5 | PPM | | 14931798 | 15024047 | 14876292 | 14872358 | |
| | | | Ave. Int. = | 14926124 | S. D. = | 70699 | | |
| Co1 | 7.00 | PPM | | 374437 | 376553 | 377215 | 379207 | |
| | | | Ave. Int. = | 378629 | S. D. = | 3879 | | |
| Co1 | 10.5 | PPM | | 25334744 | 2567195 | 25282655 | 25275709 | |
| | | | Ave. Int. = | 25345127 | S. D. = | 137283 | | |
| | | | Ave. Int. = | 25345127 | S. D. = | 137283 | | |
| Cr1 | 7.00 | PPM | | 128703 | S. D. = | 1941 | | |
| | | | Ave. Int. = | 31777 | S. D. = | 30823 | | |
| | | | Ave. Int. = | 31608 | S. D. = | 691 | | |
| | | | Ave. Int. = | 3414153 | S. D. = | 2405961 | | |
| | | | Ave. Int. = | 241607 | S. D. = | 188 | | |
| | | | Ave. Int. = | 46284 | S. D. = | 875 | | |
| | | | Ave. Int. = | 184262 | S. D. = | 10397869 | | |
| | | | Ave. Int. = | 18528518 | S. D. = | 105076 | | |
| Cr1 | 10.5 | PPM | | 7566431 | 7529145 | 7545819 | 7539495 | |
| | | | Ave. Int. = | 7565245 | S. D. = | 41943 | | |
| Cr1 | 10.5 | PPM | | 7633866 | 7616604 | 7646748 | 7608444 | |
| | | | Ave. Int. = | 7641966 | S. D. = | 32152 | | |
| Zn1 | 10.5 | PPM | | 3331651 | 3370545 | 3332467 | 3317219 | |
| | | | Ave. Int. = | 3337971 | S. D. = | 22818 | | |
| Cr1 | 10.5 | PPM | | 978703 | 981365 | 970989 | 967205 | |
| | | | Ave. Int. = | 974566 | S. D. = | 6591 | | |
| As2 | 70.0 | PPM | | 832764 | 837721 | 8323364 | 8327088 | |
| | | | Ave. Int. = | 831734 | S. D. = | 4623 | | |
| As3 | 70.0 | PPM | | 7299133 | 7271114 | 7271114 | 7271114 | |
| | | | Ave. Int. = | 72424659 | S. D. = | 40118 | | |
| As4 | 70.0 | PPM | | 3341 | 334056 | 3354718 | 3340406 | |
| | | | Ave. Int. = | 3355301 | S. D. = | 20768 | | |
| As1 | 70.0 | PPM | | 65127 | 66479 | 65481 | 65481 | |
| | | | Ave. Int. = | 65709 | S. D. = | 762 | | |
| As1 | 7.00 | PPM | | 510658 | 510657 | 5081 | 510658 | |
| | | | Ave. Int. = | 510658 | S. D. = | 343 | | |

J.E.C. Laboratories
09:56:18 17 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 13

| Line | Cone. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|-----------------------|-------|-------|--------|----------------------|------------------|-----------|----------|-------------|
| ** Standard: 4 Rep: 2 | | | | | | | | |
| As2 | 7.00 | PPM | | 102756 | 101570 | 102012 | 102892 | |
| Ba2 | 7.00 | PPM | | Ave. Int. = 26514103 | 26229236 | 26486300 | 26348763 | S. D. = 626 |
| Ca1 | 7.00 | PPM | | Ave. Int. = 26394601 | S. D. = 57203815 | 56585796 | 57245875 | 131828 |
| Cl1 | 7.00 | PPM | | Ave. Int. = 56981210 | S. D. = 11274757 | 11123552 | 11260581 | 307721 |
| Cr1 | 7.00 | PPM | | Ave. Int. = 11212163 | S. D. = 1266195 | 1245846 | 1255579 | 69800 |
| Co4 | 10.5 | PPM | | Ave. Int. = 14846936 | S. D. = 14648752 | 142738072 | 14742828 | |
| Fe1 | 7.00 | PPM | | Ave. Int. = 14757897 | S. D. = 377597 | 367053 | 375965 | 84271 |
| Ge1 | 10.5 | PPM | | Ave. Int. = 379254 | S. D. = | 372439 | 36656 | |
| Si1 | 7.00 | PPM | | Ave. Int. = 25097487 | S. D. = 252312 | 256119 | 26067 | 134519 |
| | | | | Ave. Int. = 126942 | S. D. = | 126829 | 1660 | |
| Na3 | 10.5 | PPM | | Ave. Int. = 238950 | S. D. = 238950 | 238950 | 238950 | 238950 |
| Al1 | 7.00 | PPM | | Ave. Int. = 46349 | S. D. = 46349 | 46557 | 46541 | |
| V 2 | 10.5 | PPM | | Ave. Int. = 18407227 | S. D. = 18190041 | 18335609 | 18318515 | 259 |
| Al1 | 10.5 | PPM | | Ave. Int. = 18327093 | S. D. = 18327093 | 18327093 | 18327093 | 99277 |
| Cr2 | 10.5 | PPM | | Ave. Int. = 7437806 | S. D. = 7485946 | 7503917 | 7475629 | |
| Co1 | 10.5 | PPM | | Ave. Int. = 7510524 | S. D. = 7566420 | 7584140 | 7510524 | 38236 |
| Ge1 | 10.5 | PPM | | Ave. Int. = 7566420 | S. D. = 7566420 | 7566420 | 7566420 | 7566420 |
| | | | | Ave. Int. = 3310468 | S. D. = 3310468 | 3310468 | 3310468 | 3310468 |
| As2 | 7.00 | PPM | | Ave. Int. = 8271783 | S. D. = 964202 | 964202 | 964202 | 964202 |
| | | | | Ave. Int. = 8271783 | S. D. = 817056 | 817056 | 817056 | 817056 |
| | | | | Ave. Int. = 8271783 | S. D. = 82478 | 82478 | 82478 | 82478 |
| | | | | Ave. Int. = 8271783 | S. D. = 710 | 710 | 710 | 710 |
| | | | | Ave. Int. = 8271783 | S. D. = 710 | 710 | 710 | 710 |

J.E.C. Laboratories
1255723, 17 Oct 1996

Line Conc. Units

Follows
Protocol: principle

SD/RSD

1

2

3

4

5

* Standard dilution 3

Flux

600 e.s.u.

100

PM

Ave.

Int. = 57680284 S. m. =

431278

00

AM

Int. = 31324940 S. m. =

11433834

40

PM

Int. = 31324942 S. m. =

65593

80

AM

Int. = 1266819 S. m. =

1280531

Ave. Int. = 57680284 S. m. =

1266819 S. m. =

Ave. Int. = 31324940 S. m. =

1266819 S. m. =

Ave. Int. = 31324942 S. m. =

1266819 S. m. =

Ave. Int. = 1266819 S. m. =

1266819 S. m. =

E.C. Laboratories
0:00:25 17 Oct 1996

Folder: pripli98
Protocol: prioropl

Page 15

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|----------------------|-------|-------|-------------|-----------|-----------|-----------|-----------|---|
| * Standard: 5 Rep: 1 | | | | | | | | |
| As2 | 10.0 | PPM | | 146775 | 151367 | 148921 | 146521 | |
| | | | Ave. Int. = | 148396 | S. D. = | 2254 | | |
| Br2 | 10.0 | PPM | | 37437449 | 38476551 | 37861119 | 37199938 | |
| | | | Ave. Int. = | 37743764 | S. D. = | 559858 | | |
| Ba1 | 10.0 | PPM | | 80569299 | 82758140 | 81349801 | 79979735 | |
| | | | Ave. Int. = | 81164244 | S. D. = | 1201660 | | |
| Cu1 | 10.0 | PPM | | 15748772 | 16161245 | 15916649 | 15677953 | |
| | | | Ave. Int. = | 15876155 | S. D. = | 214807 | | |
| Cr1 | 10.0 | PPM | | 1775858 | 1817633 | 1793554 | 1762026 | |
| | | | Ave. Int. = | 1787268 | S. D. = | 24006 | | |
| Cr4 | 15.0 | PPM | | 20917943 | 21435844 | 21112060 | 20784405 | |
| | | | Ave. Int. = | 21062563 | S. D. = | 282887 | | |
| Co1 | 10.0 | PPM | | 524565 | 545829 | 535405 | 525627 | |
| | | | Ave. Int. = | 532857 | S. D. = | 9930 | | |
| Mn1 | 15.0 | PPM | | 35570345 | 36523854 | 35958134 | 35395076 | |
| | | | Ave. Int. = | 35861852 | S. D. = | 500124 | | |
| Si1 | 10.0 | PPM | | 177722 | 186508 | 182714 | 178602 | |
| | | | Ave. Int. = | 181387 | S. D. = | 4049 | | |
| Tl1 | 10.0 | PPM | | 43485 | 44111 | 44031 | 42445 | |
| | | | Ave. Int. = | 43518 | S. D. = | 768 | | |
| Ni3 | 15.0 | PPM | | 3382276 | 3478325 | 3418069 | 3360738 | |
| | | | Ave. Int. = | 3409852 | S. D. = | 51409 | | |
| Sn1 | 10.0 | PPM | | 65121 | 68513 | 67519 | 65945 | |
| | | | Ave. Int. = | 66775 | S. D. = | 1527 | | |
| V2 | 15.0 | PPM | | 25946250 | 26698165 | 26270054 | 25803843 | |
| | | | Ave. Int. = | 26179578 | S. D. = | 396962 | | |
| Ca1 | 15.0 | PPM | | 10601115 | 10875984 | 10699611 | 10519815 | |
| | | | Ave. Int. = | 10674131 | S. D. = | 153339 | | |
| Fe2 | 15.0 | PPM | | 10714928 | 11009201 | 10846238 | 10678446 | |
| | | | Ave. Int. = | 10812203 | S. D. = | 149800 | | |
| Zn1 | 15.0 | PPM | | 4702630 | 4808538 | 4744078 | 4659647 | |
| | | | Ave. Int. = | 4728723 | S. D. = | 63400 | | |
| Al3 | 15.0 | PPM | | 1373483 | 1419515 | 1390333 | 1370419 | |
| | | | Ave. Int. = | 1388438 | S. D. = | 22492 | | |
| Na2 | 100. | PPM | | 1176708 | 1211168 | 1191658 | 1170682 | |
| | | | Ave. Int. = | 1187554 | S. D. = | 18044 | | |
| Cl3 | 100. | PPM | | 101341374 | 104056354 | 102426379 | 100615604 | |
| | | | Ave. Int. = | 102109928 | S. D. = | 1495816 | | |
| Mg4 | 100. | PPM | | 4738834 | 4862178 | 4791730 | 4714346 | |
| | | | Ave. Int. = | 4776772 | S. D. = | 65458 | | |
| K1 | 100. | PPM | | 93661 | 95435 | 93471 | 92221 | |
| | | | Ave. Int. = | 93697 | S. D. = | 1323 | | |
| Mo1 | 10.0 | PPM | | 722390 | 748542 | 733614 | 721198 | |
| | | | Ave. Int. = | 731436 | S. D. = | 12702 | | |

S.E.C. Laboratories
10:01:22 17 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 16

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|-----------------------|-------|-------|-------------|----------|--------------------------|-----------|----------|---|
| ** Standard: 5 Rep: 2 | | | | | | | | |
| | | | | Seq: 537 | 10:01:22 17 Oct 1996 ICP | | | |
| As2 | 10.0 | PPM | | 144449 | 145607 | 144441 | 141879 | |
| | | | Ave. Int. = | 144094 | S. D. = | 1575 | | |
| Se2 | 10.0 | PPM | | 36800922 | 37007662 | 37083200 | 36132748 | |
| | | | Ave. Int. = | 36756133 | S. D. = | 432379 | | |
| Ba1 | 10.0 | PPM | | 79013660 | 79512832 | 79471105 | 77365443 | |
| | | | Ave. Int. = | 78840760 | S. D. = | 1009203 | | |
| Dd1 | 10.0 | PPM | | 15542869 | 15571503 | 15670103 | 15242824 | |
| | | | Ave. Int. = | 15506825 | S. D. = | 184245 | | |
| Co1 | 10.0 | PPM | | 1745690 | 1748260 | 1753010 | 1706718 | |
| | | | Ave. Int. = | 1738420 | S. D. = | 21351 | | |
| Cr4 | 15.0 | PPM | | 20584086 | 20650479 | 20716428 | 20137383 | |
| | | | Ave. Int. = | 20522094 | S. D. = | 262103 | | |
| pb1 | 10.0 | PPM | | 520629 | 528269 | 524877 | 514621 | |
| | | | Ave. Int. = | 522099 | S. D. = | 5884 | | |
| Mn1 | 15.0 | PPM | | 35023260 | 35181834 | 35270301 | 34383454 | |
| | | | Ave. Int. = | 34964712 | S. D. = | 400756 | | |
| b1 | 10.0 | PPM | | 177840 | 179506 | 177870 | 175834 | |
| | | | Ave. Int. = | 177763 | S. D. = | 1503 | | |
| Tl1 | 10.0 | PPM | | 42813 | 42951 | 43893 | 42511 | |
| | | | Ave. Int. = | 43042 | S. D. = | 596 | | |
| Ni3 | 15.0 | PPM | | 3331094 | 3359290 | 3359340 | 3262056 | |
| | | | Ave. Int. = | 3327945 | S. D. = | 45896 | | |
| Te1 | 10.0 | PPM | | 64839 | 64913 | 65049 | 63391 | |
| | | | Ave. Int. = | 64548 | S. D. = | 776 | | |
| V2 | 15.0 | PPM | | 25548406 | 25640135 | 25747964 | 25059201 | |
| | | | Ave. Int. = | 25498927 | S. D. = | 304284 | | |
| Uu1 | 15.0 | PPM | | 10428855 | 10471029 | 10485233 | 10201133 | |
| | | | Ave. Int. = | 10396563 | S. D. = | 132468 | | |
| Fe2 | 15.0 | PPM | | 10563649 | 10604571 | 10627941 | 10356734 | |
| | | | Ave. Int. = | 10538224 | S. D. = | 123876 | | |
| Sn1 | 15.0 | PPM | | 4635469 | 4642909 | 4653839 | 4535810 | |
| | | | Ave. Int. = | 4617007 | S. D. = | 54654 | | |
| I13 | 15.0 | PPM | | 1354501 | 1360661 | 1365077 | 1337869 | |
| | | | Ave. Int. = | 1354527 | S. D. = | 11922 | | |
| Na2 | 100. | PPM | | 1158656 | 1162754 | 1169722 | 1127985 | |
| | | | Ave. Int. = | 1154779 | S. D. = | 18438 | | |
| a3 | 100. | PPM | | 99790435 | 100126312 | 100405879 | 97755148 | |
| | | | Ave. Int. = | 99519444 | S. D. = | 1202807 | | |
| Mg4 | 100. | PPM | | 4673130 | 4699292 | 4703394 | 4581406 | |
| | | | Ave. Int. = | 4664306 | S. D. = | 56869 | | |
| Li1 | 100. | PPM | | 92787 | 92583 | 91649 | 90883 | |
| | | | Ave. Int. = | 91976 | S. D. = | 881 | | |
| Mo1 | 10.0 | PPM | | 714710 | 712598 | 720040 | 697640 | |
| | | | Ave. Int. = | 711247 | S. D. = | 9597 | | |

S.E.C. Laboratories
10:02:19 17 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 17

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|------------------------|-------|-------|-------------|-----------|-----------|-----------|-----------|---|
| *** Standard: 5 Rep: 3 | | | | | | | | |
| As2 | 10.0 | PPM | | 148329 | 148289 | 147303 | 150383 | |
| | | | Ave. Int. = | 148576 | S. D. = | 1295 | | |
| Ba2 | 10.0 | PPM | | 37269906 | 37360934 | 37173151 | 37772367 | |
| | | | Ave. Int. = | 37394090 | S. D. = | 263583 | | |
| Ba1 | 10.0 | PPM | | 79967111 | 79983673 | 79733934 | 81244088 | |
| | | | Ave. Int. = | 80232202 | S. D. = | 684160 | | |
| Cd1 | 10.0 | PPM | | 15695360 | 15764732 | 15646769 | 15882569 | |
| | | | Ave. Int. = | 15747358 | S. D. = | 102316 | | |
| Co1 | 10.0 | PPM | | 1763548 | 1772010 | 1761658 | 1791554 | |
| | | | Ave. Int. = | 1772193 | S. D. = | 13670 | | |
| Cr4 | 15.0 | PPM | | 20830827 | 20881773 | 20770439 | 21089660 | |
| | | | Ave. Int. = | 20893175 | S. D. = | 138670 | | |
| Fe1 | 10.0 | PPM | | 528749 | 528485 | 532381 | 533269 | |
| | | | Ave. Int. = | 530721 | S. D. = | 2459 | | |
| Mn1 | 15.0 | PPM | | 35432908 | 35534746 | 35396828 | 35925527 | |
| | | | Ave. Int. = | 35572502 | S. D. = | 242487 | | |
| Se1 | 10.0 | PPM | | 183620 | 178278 | 177882 | 181094 | |
| | | | Ave. Int. = | 180219 | S. D. = | 2681 | | |
| Tl1 | 10.0 | PPM | | 44349 | 43557 | 44459 | 44179 | |
| | | | Ave. Int. = | 44136 | S. D. = | 403 | | |
| Na3 | 15.0 | PPM | | 3368594 | 3377532 | 3363708 | 3422931 | |
| | | | Ave. Int. = | 3383191 | S. D. = | 27104 | | |
| Se1 | 10.0 | PPM | | 65657 | 65393 | 65865 | 66257 | |
| | | | Ave. Int. = | 65793 | S. D. = | 365 | | |
| V2 | 15.0 | PPM | | 25893852 | 25973164 | 25821595 | 26283196 | |
| | | | Ave. Int. = | 25992952 | S. D. = | 203156 | | |
| U1 | 15.0 | PPM | | 10555621 | 10573119 | 10527899 | 10713677 | |
| | | | Ave. Int. = | 10592579 | S. D. = | 82851 | | |
| Fe2 | 15.0 | PPM | | 10689528 | 10711656 | 10651217 | 10847420 | |
| | | | Ave. Int. = | 10724955 | S. D. = | 85375 | | |
| Sn1 | 15.0 | PPM | | 4681149 | 4679631 | 4655061 | 4739318 | |
| | | | Ave. Int. = | 4688790 | S. D. = | 35744 | | |
| Al3 | 15.0 | PPM | | 1377757 | 1378925 | 1371259 | 1392827 | |
| | | | Ave. Int. = | 1380192 | S. D. = | 9073 | | |
| Na2 | 100. | PPM | | 1172026 | 1172688 | 1170208 | 1192264 | |
| | | | Ave. Int. = | 1176797 | S. D. = | 10365 | | |
| Ka3 | 100. | PPM | | 100872279 | 101101733 | 100491329 | 102009070 | |
| | | | Ave. Int. = | 101118603 | S. D. = | 644816 | | |
| Mg4 | 100. | PPM | | 4721610 | 4743002 | 4711892 | 4772658 | |
| | | | Ave. Int. = | 4737291 | S. D. = | 26922 | | |
| Li1 | 100. | PPM | | 92839 | 93057 | 93375 | 94589 | |
| | | | Ave. Int. = | 93465 | S. D. = | 781 | | |
| Mo1 | 10.0 | PPM | | 723336 | 723334 | 719824 | 737238 | |
| | | | Ave. Int. = | 725933 | S. D. = | 7716 | | |

E.C. Laboratories
10:08:02 17 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 18

| Name | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|------------------------|-----------|-------------|-------------|----------|--------------------------|--------|---|---|
| *** Standard: 1 Rep: 1 | | | | Seq: 539 | 10:08:02 17 Oct 1996 ICP | | | |
| Ag1 | .000 | PPM | -1441 | -2197 | -1443 | -1285 | | |
| | | | Ave. Int. = | -1592 | S. D. = | 410 | | |
| *** Standard: 1 Rep: 2 | | | | Seq: 540 | 10:08:51 17 Oct 1996 ICP | | | |
| Ag1 | .000 | PPM | -677 | 2215 | 871 | -1657 | | |
| | | | Ave. Int. = | 188 | S. D. = | 1706 | | |
| *** Standard: 1 Rep: 3 | | | | Seq: 541 | 10:09:39 17 Oct 1996 ICP | | | |
| Ag1 | .000 | PPM | -465 | -477 | -1215 | 795 | | |
| | | | Ave. Int. = | -341 | S. D. = | 834 | | |
| ** Peak | | | | Seq: 542 | 10:11:19 17 Oct 1996 | | | |
| XAp = 215 | YAp = 141 | Intensity = | 7548068 | | | | | |
| *** Standard: 5 Rep: 1 | | | | Seq: 543 | 10:12:39 17 Oct 1996 ICP | | | |
| Ag1 | 1.00 | PPM | 270721 | 257215 | 236825 | 226957 | | |
| | | | Ave. Int. = | 247930 | S. D. = | 19739 | | |
| *** Standard: 5 Rep: 2 | | | | Seq: 544 | 10:13:27 17 Oct 1996 ICP | | | |
| Ag1 | 1.00 | PPM | 218617 | 224581 | 235097 | 241149 | | |
| | | | Ave. Int. = | 229861 | S. D. = | 10151 | | |
| *** Standard: 5 Rep: 3 | | | | Seq: 545 | 10:14:15 17 Oct 1996 ICP | | | |
| Ag1 | 1.00 | PPM | 248771 | 255611 | 245123 | 252071 | | |
| | | | Ave. Int. = | 250394 | S. D. = | 4489 | | |
| *** Standard: 6 Rep: 1 | | | | Seq: 546 | 10:17:16 17 Oct 1996 ICP | | | |
| Ag1 | 2.00 | PPM | 543690 | 554580 | 561034 | 534154 | | |
| | | | Ave. Int. = | 548365 | S. D. = | 11873 | | |
| *** Standard: 6 Rep: 2 | | | | Seq: 547 | 10:18:04 17 Oct 1996 ICP | | | |
| Ag1 | 2.00 | PPM | 514840 | 513698 | 477804 | 443827 | | |
| | | | Ave. Int. = | 487542 | S. D. = | 33839 | | |
| *** Standard: 6 Rep: 3 | | | | Seq: 548 | 10:18:52 17 Oct 1996 ICP | | | |
| Ag1 | 2.00 | PPM | 466842 | 495836 | 482780 | 441227 | | |
| | | | Ave. Int. = | 471671 | S. D. = | 23505 | | |

E.C. Laboratories
10:21:53 17 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 19

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|---------------------------|-------|-------|--------|------|-------|--------------------------|-------|---|
| *** Update Intercept | | | | | | | | |
| | | | Corr. | Seq: | 549 | 10:21:53 17 Oct 1996 ICP | | |
| Ag1 | .000 | PPM | -.025 | Avg. | .025 | .006 | .025 | |
| As2 | .000 | PPM | .028 | | -.028 | .038 | -.028 | |
| Br2 | .000 | PPM | .034 | | -.034 | .000 | -.034 | |
| Cd1 | .000 | PPM | .038 | | -.038 | .001 | -.038 | |
| Cl1 | .000 | PPM | .068 | | -.068 | .001 | -.068 | |
| Cr4 | .000 | PPM | .057 | | -.057 | .005 | -.057 | |
| Pb1 | .000 | PPM | .078 | | -.078 | .003 | -.078 | |
| Mn1 | .000 | PPM | .056 | | -.056 | .019 | -.056 | |
| Si1 | .000 | PPM | .061 | | -.061 | .001 | -.061 | |
| Tl1 | .000 | PPM | -.013 | | .013 | .038 | .013 | |
| Ni3 | .000 | PPM | .041 | | -.041 | .061 | -.041 | |
| S11 | .000 | PPM | .092 | | -.092 | .005 | -.092 | |
| V2 | .000 | PPM | -.008 | | .008 | .036 | .008 | |
| Cu1 | .000 | PPM | .048 | | -.048 | .002 | -.048 | |
| F2 | .000 | PPM | .057 | | -.057 | .002 | -.057 | |
| Zn1 | .000 | PPM | .073 | | -.073 | .005 | -.073 | |
| Al3 | .000 | PPM | .064 | | -.064 | .002 | -.064 | |
| Na2 | .000 | PPM | .032 | | -.032 | .015 | -.032 | |
| Cl3 | .000 | PPM | .234 | | -.234 | .114 | -.234 | |
| Mg4 | .000 | PPM | .424 | | -.424 | .001 | -.424 | |
| K1 | .000 | PPM | .130 | | -.130 | .033 | -.130 | |
| Mg1 | .000 | PPM | .067 | | -.067 | .266 | -.067 | |
| | | | .029 | | -.029 | .009 | -.029 | |
| *** Check Standard: 2 Ck2 | | | | | | | | |
| Line | Flag | %Rcv. | Found | True | Units | SD/RSD | | |
| Ag1 | | 101. | 1.01 | 1.00 | PPM | .090 | | |
| Ba1 | | 107. | 1.07 | 1.00 | PPM | .014 | | |
| Al3 | | 107. | 1.07 | 1.00 | PPM | .031 | | |

E.C. Laboratories
10:30:56 17 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 20

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|---------------------------|---------|-----------|-------------|------|---------|----------------------|---|---|
| *** Check Standard: 3 Ck3 | | | | | | | | |
| Line | Flag | %Rcv. | Found | True | Units | SD/RSD | | |
| As2 | H | 111. | 1.11 | 1.00 | PPM | .060 | | |
| Be2 | | 105. | 1.05 | 1.00 | PPM | .014 | | |
| Cu1 | | 107. | 1.07 | 1.00 | PPM | .010 | | |
| Co1 | H | 111. | 1.11 | 1.00 | PPM | .020 | | |
| Cr4 | H | 112. | 1.12 | 1.00 | PPM | .015 | | |
| Pb1 | | 108. | 1.08 | 1.00 | PPM | .042 | | |
| Mn1 | | 107. | 1.07 | 1.00 | PPM | .012 | | |
| Sb1 | | 98.6 | .986 | 1.00 | PPM | .018 | | |
| Tl1 | | 98.3 | .983 | 1.00 | PPM | .124 | | |
| Ni3 | H | 112. | 1.12 | 1.00 | PPM | .009 | | |
| Se1 | H | 113. | 1.13 | 1.00 | PPM | .055 | | |
| V2 | | 109. | 1.09 | 1.00 | PPM | .013 | | |
| Ca1 | | 106. | 1.06 | 1.00 | PPM | .011 | | |
| Fe2 | H | 111. | 1.11 | 1.00 | PPM | .008 | | |
| Zn1 | | 108. | 1.08 | 1.00 | PPM | .014 | | |
| Mo1 | H | 113. | 1.13 | 1.00 | PPM | .034 | | |
| *** Peak | | | | | | | | |
| | | | | Seq: | 552 | 10:33:16 17 Oct 1996 | | |
| X | p = 213 | YAp = 140 | Intensity = | | 7576852 | | | |
| *** Check Standard: 4 Ck4 | | | | | | | | |
| Line | Flag | %Rcv. | Found | True | Units | SD/RSD | | |
| Mg2 | | 100. | 50.2 | 50.0 | PPM | .630 | | |
| Ca3 | | 103. | 51.7 | 50.0 | PPM | .739 | | |
| Mn4 | | 103. | 10.3 | 10.0 | PPM | .163 | | |
| Kr1 | | 106. | 10.6 | 10.0 | PPM | .460 | | |
| *** Check Standard: 3 Ck3 | | | | | | | | |
| Line | Flag | %Rcv. | Found | True | Units | SD/RSD | | |
| As2 | | 108. | 1.08 | 1.00 | PPM | .096 | | |
| Be2 | | 107. | 1.07 | 1.00 | PPM | .017 | | |
| Cu1 | H | 113. | 1.13 | 1.00 | PPM | .023 | | |
| Co1 | | 107. | 1.07 | 1.00 | PPM | .027 | | |
| Cr4 | | 105. | 1.05 | 1.00 | PPM | .020 | | |
| Pb1 | H | 111. | 1.11 | 1.00 | PPM | .019 | | |
| Mn1 | | 108. | 1.08 | 1.00 | PPM | .018 | | |
| Sb1 | | 102. | 1.02 | 1.00 | PPM | .049 | | |
| Tl1 | H | 110. | 1.10 | 1.00 | PPM | .035 | | |
| Ni3 | | 102. | 1.02 | 1.00 | PPM | .014 | | |
| Se1 | | 107. | 1.07 | 1.00 | PPM | .041 | | |
| V2 | | 107. | 1.07 | 1.00 | PPM | .020 | | |
| Ca1 | H | 110. | 1.10 | 1.00 | PPM | .016 | | |
| Fe2 | | 105. | 1.05 | 1.00 | PPM | .022 | | |
| Zn1 | | 109. | 1.09 | 1.00 | PPM | .015 | | |
| Mo1 | | 107. | 1.07 | 1.00 | PPM | .036 | | |

E.C. Laboratories
10:46:03 17 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 21

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|-------|-------|--------------|---------|------|--------|--------------------------|---|---|
| ----- | | | | | | | | |
| * | * | Update Slope | | Seq: | 555 | 10:46:03 17 Oct 1996 ICP | | |
| | | | % Corr. | Avg. | SD/RSD | Old | | |
| As2 | 3.00 | PPM | -4.42 | 3.14 | .060 | 3.14 | | |
| Be2 | 3.00 | PPM | -3.37 | 3.10 | .010 | 3.10 | | |
| B_1 | 3.00 | PPM | -2.95 | 3.09 | .011 | 3.09 | | |
| Ca1 | 3.00 | PPM | -5.23 | 3.17 | .009 | 3.17 | | |
| Co1 | 3.00 | PPM | -4.39 | 3.14 | .011 | 3.14 | | |
| C_4 | 4.50 | PPM | -4.71 | 4.72 | .015 | 4.72 | | |
| P_1 | 3.00 | PPM | -3.11 | 3.10 | .024 | 3.10 | | |
| Mn1 | 4.50 | PPM | -3.48 | 4.66 | .016 | 4.66 | | |
| Sb1 | 3.00 | PPM | -4.87 | 3.15 | .066 | 3.15 | | |
| Tl_1 | 3.00 | PPM | -4.19 | 3.13 | .109 | 3.13 | | |
| Ni3 | 4.50 | PPM | -3.92 | 4.68 | .026 | 4.68 | | |
| Se1 | 3.00 | PPM | -7.14 | 3.23 | .096 | 3.23 | | |
| V_2 | 4.50 | PPM | -3.39 | 4.66 | .018 | 4.66 | | |
| Ca_1 | 4.50 | PPM | -2.86 | 4.63 | .014 | 4.63 | | |
| Fe2 | 4.50 | PPM | -4.39 | 4.71 | .019 | 4.71 | | |
| Zn_1 | 4.50 | PPM | -4.02 | 4.69 | .010 | 4.69 | | |
| A_3 | 4.50 | PPM | -2.49 | 4.61 | .035 | 4.61 | | |
| Na2 | 30.0 | PPM | -2.11 | 30.6 | .245 | 30.6 | | |
| Ca3 | 30.0 | PPM | -2.18 | 30.7 | .106 | 30.7 | | |
| Mg_4 | 30.0 | PPM | -3.63 | 31.1 | .109 | 31.1 | | |
| K_1 | 30.0 | PPM | -4.17 | 31.3 | .423 | 31.3 | | |
| Mo1 | 3.00 | PPM | -6.82 | 3.22 | .016 | 3.22 | | |
| * | * | Update Slope | | Seq: | 556 | 10:49:09 17 Oct 1996 ICP | | |
| | | | % Corr. | Avg. | SD/RSD | Old | | |
| A_11 | 1.00 | PPM | .420 | .996 | .071 | .996 | | |

* * Check Standard: 2 Ck2 Seq: 557 10:54:05 17 Oct 1996 ICP
 Line Flag %Rcv. Found True Units SD/RSD
 401 102. 1.02 1.00 PPM .041
 31 103. 1.03 1.00 PPM .010
 413 103. 1.03 1.00 PPM .015
 * * Peak Seq: 558 10:55:45 17 Oct 1996
 XAp = 212 YAp = 140 Intensity = 7592380

E.C. Laboratories
10:57:05 17 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 22

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|---------------------------|-------|-------|--------|------|-------|----------|--------------------------|--------|
| *** Check Standard: 3 Ck3 | | | | | | | | |
| As2 | 99.1 | .991 | | True | Units | Seq: 559 | 10:57:05 17 Oct 1996 ICP | SD/RSD |
| Br2 | 99.2 | .992 | | 1.00 | PPM | | | .032 |
| Co1 | 100. | 1.00 | | 1.00 | PPM | | | .003 |
| Co1 | 101. | 1.01 | | 1.00 | PPM | | | .011 |
| Cr4 | 98.6 | .986 | | 1.00 | PPM | | | .015 |
| Cr1 | 99.5 | .995 | | 1.00 | PPM | | | .003 |
| Mn1 | 97.7 | .977 | | 1.00 | PPM | | | .018 |
| Sb1 | 100. | 1.00 | | 1.00 | PPM | | | .005 |
| Tl1 | 102. | 1.02 | | 1.00 | PPM | | | .029 |
| Ni3 | 96.1 | .961 | | 1.00 | PPM | | | .057 |
| Se1 | 98.1 | .981 | | 1.00 | PPM | | | .010 |
| V2 | 97.9 | .979 | | 1.00 | PPM | | | .093 |
| Cl1 | 101. | 1.01 | | 1.00 | PPM | | | .006 |
| Fe2 | 98.2 | .982 | | 1.00 | PPM | | | .005 |
| Zn1 | 96.9 | .969 | | 1.00 | PPM | | | .008 |
| Mg1 | 97.1 | .971 | | 1.00 | PPM | | | .008 |
| *** Check Standard: 4 Ck4 | | | | | | | | |
| Na2 | 98.9 | 49.4 | | True | Units | Seq: 560 | 11:00:08 17 Oct 1996 ICP | SD/RSD |
| Ca3 | 101. | 50.6 | | 50.0 | PPM | | | .330 |
| Mg4 | 100. | 10.0 | | 50.0 | PPM | | | .260 |
| K1 | 102. | 10.2 | | 10.0 | PPM | | | .074 |
| | | | | 10.0 | PPM | | | .559 |

E.C. Laboratories
11:02:40 17 Oct 1996

Folder: pripl198
Protocol: priorplt

Page 23

| Line | Conc. | Units | SD/RSD | 1 | 2 | 3 | 4 | 5 |
|--------------------|-------|-------|--------|-------|--------|-------|---|---|
| * Update Intercept | | | | | | | | |
| | | | Corr. | Avg. | SD/RSD | Old | | |
| Ag1 | .000 | PPM | -.020 | .020 | .006 | .020 | | |
| Zn2 | .000 | PPM | .018 | -.018 | .019 | -.018 | | |
| Bi2 | .000 | PPM | .034 | -.033 | .000 | -.033 | | |
| Ba1 | .000 | PPM | .003 | -.003 | .015 | -.003 | | |
| Cd1 | .000 | PPM | .070 | -.066 | .001 | -.066 | | |
| Y1 | .000 | PPM | .061 | -.059 | .005 | -.059 | | |
| Ca4 | .000 | PPM | .079 | -.075 | .002 | -.075 | | |
| Pb1 | .000 | PPM | .064 | -.062 | .013 | -.062 | | |
| V1 | .000 | PPM | .061 | -.058 | .001 | -.058 | | |
| Li1 | .000 | PPM | .048 | -.046 | .072 | -.046 | | |
| Tl1 | .000 | PPM | .078 | -.075 | .043 | -.075 | | |
| Ni3 | .000 | PPM | .089 | -.086 | .004 | -.086 | | |
| S1 | .000 | PPM | -.019 | .018 | .046 | .018 | | |
| V2 | .000 | PPM | .055 | -.053 | .002 | -.053 | | |
| Cu1 | .000 | PPM | .063 | -.061 | .003 | -.061 | | |
| F2 | .000 | PPM | .075 | -.071 | .003 | -.071 | | |
| Zn1 | .000 | PPM | .062 | -.060 | .003 | -.060 | | |
| Al3 | .000 | PPM | -.020 | .019 | .052 | .019 | | |
| Ni2 | .000 | PPM | .194 | -.190 | .180 | -.190 | | |
| Ca3 | .000 | PPM | .387 | -.379 | .015 | -.379 | | |
| Mg4 | .000 | PPM | .102 | -.098 | .024 | -.098 | | |
| K1 | .000 | PPM | .014 | -.013 | .386 | -.013 | | |
| Mn1 | .000 | PPM | .003 | -.003 | .017 | -.003 | | |

Section VI
General Chemistry Raw Data & QC Reports

Quality Control Data

| | | | |
|---------------------------|---|--------------------------------|------------------|
| Client: | <u>Ecology and Environment</u> | | |
| Project : | <u>S05-9610-805</u> | Date Sampled | <u>10/11/96*</u> |
| Project No.: | <u>KJ5100</u> | Date Received in Lab | <u>10/15/96</u> |
| Sample Types: | <u>Vat & Drum Samples</u> | Date Analysis Completed | <u>10/17/96</u> |
| Station No.: | <u>Drum 01, 02, 04, 10, 11 & Vat 08.</u> | | |
| B.E.C. Sample No.: | <u>96T14933, 96T14934, 96T14936, 96T14940, 96T14942 & 96T14943.</u> | | |

| Test | Batch No. | Method Blank | Method Spike | Duplicate | | RPD % | Matrix Spike | units |
|------------------|---|-----------------|---|----------------|----------------|----------|-----------------|-------|
| | | | | R ₁ | R ₂ | | | |
| Total Chloride | 457-45 | LT 1 mg/l | 101.48% | 19.2 | 19.2 | 0.0% | 100.32% | mg/l |
| Total Sulfate | 457-48 | LT 2 mg/l | 96.7% | 1720 | 1740 | 1.16% | 95.52% | mg/l |
| Total Cyanide | 446-44 | LT 0.1 | 82.16% | 2.79 | 2.82 | 1.07% | 100.1% | mg/kg |
| Reactive Cyanide | 446-47 | LT 0.1 | 22.26% | LT 0.1 | LT 0.1 | NA | 0.6%** | mg/kg |
| Flash Point | | NA | Calibrated on 10/7/96 using p-xylene 79.0°F (Limits are 81 ± 1.5°F). | | | | | |
| pH | Meter calibrated on 10/15/96 buffer values 3.99, 7.00 & 10.01 | | (Limits ±0.05 S.U. from true values). | | | | | S.U. |

pH meter was standardized on the day of analysis.

* Date sampled per chain-of-custody

** Low recovery for Reactive Cyanide is typical for this procedure.

Approved By: J. F. Klein
 Date: 10/25/96

Biological & Environmental Control Laboratories Inc.

Chemistry Department

RAW DATA RESULT LOG

Analytical Test:

Reactive CN

Batch No: 446 - 47

Batch Start Date: 10/15/96

Matrix:

Batch End Date:

| | SAMPLE # | weight(g) | m1 dist. filtr | Abu + 578nm | | DATE ANALYZED | INITIALS |
|-------|-----------------|------------------|-----------------------|--------------------|--|----------------------|-----------------|
| MBLK | QC sand | 10.12 | 25 | 0.001 | | 10/15/96 | CSY |
| MSPK | QC sand | 10.48 | 25 | 0.1844 0.652 | | | |
| 1 | 96T14940 | 10.57 | 25 | 0.027 0.068 | | | |
| DUP | 96T14940 | 10.57 | 25 | 0.011 | | | |
| SPIKE | 96T14940 | 10.57 | 25 | 0.027 0.068 | | J | |
| 2 | 96T14943 | 10.07 | 25 | 0.554 | | 10/15/96 | CSY |
| 3 | | | | | | | |
| 4 | | | | | | | |
| 5 | | | | | | | |
| 6 | | | | | | | |
| 7 | | | | | | | |
| 8 | | | | | | | |
| 9 | | | | | | | |
| 10 | | | | | | | |

COMMENTS:

$S_{\text{spike}} = 2.5 \text{ mg} / 2.8 \text{ g}$
100%

QA/QC:

METHOD BLANK = LT 0.1 mg/kg

| SAMPLE | AMT ADDED | AMT RECOV. | % RECOV. |
|---------------|------------------|-------------------|-----------------|
| METH SPK. | 200 mg | 45.2 mg | 22.6 |
| R1 (RESULT) | | 0.0565 mg | |
| R2 (RESULT) | | 0.0758 mg | 24 %RPD |
| MATRIX SPK | 200 mg | 12 mg | 0.6 |

Approved by QA/QC: S.F.K 10/20/96

+ Sample Dip at lot 100 mg

Page

Biological & Environmental Control Laboratories Inc.

Chemistry Department

RAW DATA RESULT LOG

Batch No: 4416-414

Batch Start Date: 10-7-96

Matrix: Solids

Batch End Date:

| | SAMPLE # | wt. Sample | m/Dist | ABS @ 578nm | DATE ANALYZED | INITIALS |
|-------|----------|----------------|------------------------------|--------------------|---------------|----------|
| MBLK | QCSand | 10.10 | 25 | 0.005 0.007 cmv | 10-7-96 | cmv |
| MSPK | QCSand | 10.06 | 25 | 0.248 | 10-7-96 | cmv |
| 1 | 96T14533 | 10.02 | 25 | 0.429 | 10-8-96 | |
| DUP | 96T14533 | 10.07 | 25 | 0.435 | | |
| SPIKE | 96T14533 | 10.05 | 25 | 0.746 | | |
| 2 | 96T14534 | 10.15 | 15 | 0.447 | ↓ | ↓ |
| 3 | 96T14537 | 10.08 | 25 | 0.045 | 10-8-96 | cmv |
| 4 | 96T14689 | 10.22 10.40 | 25 | 0.062 | 10-9-96 | cmv |
| 5 | 96T14781 | 10.19 | 25 | 0.007 | 10-14-96 | cmv |
| 6 | 96T14782 | 10.06 | 25 | 0.048 | 10-14-96 | cmv |
| 7 | 96T14940 | 10.27 | 125 ± .2ml used = 0.003ml | 0.323 | 10-16-96 | |
| 8 | 96T14943 | 10.47 | 150 ± .1ml used = 0.002ml | .693 | ↓ | ↓ |
| 9 | 96T14984 | 10.34 | 25 | 0.307 | 10-16-96 | cmv |
| 10 | | | | | | |

COMMENTS:

Spk IOC - 259 2ml 10mg/L
 rSTD (6ug) = 0.354 10-7-96
 " " = 0.345 10-8-96
 " " = 0.351 10-9-96
 " " = 0.357 10-14-96
 " " = 0.321 10-16-96

QA/QC:

| SAMPLE | AMT ADDED | AMT RECOV | % RECOV |
|------------------------------------|-----------|-------------|-----------|
| METH SPK. | 20ug | 16.4ug | 82.16 |
| R1 (RESULT) | | 2.793 mg/kg | |
| R2 (RESULT) | | 2.817 mg/kg | .86 % RPD |
| MATRIX SPK | 20ug | 20.0ug | 100.14 |
| Approved by QA/QC: J.F.B. 10/20/96 | | | |

Biological & Environmental Control Laboratories Inc.

Chemistry Department

Analytical Test:

²
SO₄

RAW DATA RESULT LOG

Batch No: 457-48

Batch Start Date: 10/15/96

Matrix: Liquid

Batch End Date:

| | SAMPLE # | ml Sample | NTU | | | DATE ANALYZED | INITIALS |
|-------|----------------------|---|-------|--|--|---------------|----------|
| MBLK | D.I H ₂ O | 100 | 0.42 | | | 10/15/96 | AAT |
| MSPK | D.I H ₂ O | 100 | 94.5 | | | | |
| 1 | 96T14760 | 1.0 | 83.7 | | | | |
| DUP | 96T14760 | 1.0 | 84.8 | | | | |
| SPIKE | 96T14760 | 1.0 | 179.4 | | | V | V |
| 2 | 96T14933 | 1/100 > 1/10,000 used 20 ml = 0.002 ml | 149.4 | | | 10/15/96 | AAT |
| R2 | 96T14933 | 1/1000 used 2 ml = 0.002 ml | 152.3 | | | 10/16/96 | AAT |
| 4 | 96T15120 | 2.0 | 173.2 | | | 10/21/96 | AAT |
| 5 | | | | | | | |
| 6 | | | | | | | |
| 7 | | | | | | | |
| 8 | | | | | | | |
| 9 | | | | | | | |
| 10 | | | | | | | |

COMMENTS:

Spk 10ml ± 528

200 ppm

✓ curve 98.2

✓ 98.9

✓ 96.4

QA/QC:

METHOD BLANK = LT 2.0 mg/L

| SAMPLE | AMT ADDED | AMT RECOV. | % RECOV. |
|-------------|-----------|------------|------------|
| METH SPK. | 2000 µg | 1934 µg | 96.7 |
| R1 (RESULT) | | 1717.74 µg | |
| R2 (RESULT) | | 1739.82 µg | 1.28 % RPD |
| MATRIX SPK | 2000 µg | 1910 µg | 95.52 |

Approved by QA/QC: JFB 10/20/96

Page

0048

Biological & Environmental Control Laboratories Inc.
Chemistry Department

Analytical Test:

CHLORINE

RAW DATA RESULT LOG

Batch No: 457-45

Batch Start Date: 10-14-96

Matrix:

Batch End Date:

| | SAMPLE # | ml Sample | ml Hg(No ₃) ₂ | | DATE ANALYZED | INITIALS |
|-------|---------------------|-----------|--------------------------------------|--|---------------|----------|
| MBLK | DI H ₂ O | .50 | 0.05 | | 10-14-96 | GJB |
| MSPK | ↓ | .50 | 8.85 | | | ↓ |
| 1 | 96T14760 | .50 | 1.80 | | | ↓ |
| DUP | 96T14760 | .50 | 1.80 | | | ↓ |
| SPIKE | 96T14760 | .50 | 10.45 14.45 14.46 | | | ↓ |
| 2 | 96T14934 | 0.1 ml | 70.1 ml | | 10-17-96 | GJB |
| 3 | 96T14936 | .0075 ml | 2.4 ml | | ↓ | ↓ |
| 4 | 96T15120 | .50 | 1.7 | | 10/26/96 | |
| 5 | | | | | | |
| 6 | | | | | | |
| 7 | | | | | | |
| 8 | | | | | | |
| 9 | | | | | | |
| 10 | | | | | | |

COMMENTS:

SD, Hg = 5.0 mg & 531

No. molality Hg(No₃)₂ = 17.3 ml

No. molar Hg(No₃)₂ = 16.2 ml

N Hg(No₃)₂ = 17.7 10/26/96

QA/QC:

METHOD BLANK = LTL mg/l

| SAMPLE | AMT ADDED | AMT RECOV. | % RECOV. |
|-------------|-----------|------------|-----------|
| METH SPK. | 88.6 µg | 90.42 µg | 101.48% |
| R1 (RESULT) | | 19.2 mg/l | |
| R2 (RESULT) | | 19.2 mg/l | 0.00% RPD |
| MATRIX SPK | 88.6 µg | 106.86 µg | 110.32% |

Approved by QA/QC: J.F.B. 10/26/96

Month: OCTOBER 1996

ORION SA720 pH METER

| | Date Checked | pH 4.00 | pH 7.00 | pH 10.00 | Temp | Slope | Initials |
|----|--------------|---------|---------|----------|--------|--------|----------|
| 1 | 10/1/96 | 4.00 | 7.00 | 9.97 | 22.6 | 96.5 | FG |
| 2 | 10/2/96 | 4.01 | 7.00 | 9.99 | 22.7 | 96.8 | JEB |
| 3 | 10/3/96 | 4.01 | 7.00 | 9.96 | 22.2°C | 98.1 | JEB/HAT |
| 4 | | | | | | | |
| 5 | | | | | | | |
| 6 | | | | | | | |
| 7 | 10/7/96 | 4.01 | 6.99 | 10.02 | 21.8°C | 98.5 | JEB |
| 8 | 10/8/96 | 4.01 | 7.00 | 10.01 | 22.4°C | 98.5 | JEB |
| 9 | 10/9/96 | 4.01 | 7.00 | 9.99 | 22.4°C | 98.2 | JEB |
| 10 | 10/10/96 | 4.00 | 7.01 | 10.00 | 22.4°C | 98.9 | JEB |
| 11 | 10/11/96 | 4.00 | 7.00 | 10.00 | 22.1 | 97.7 | CY |
| 12 | | | | | | | |
| 13 | | | | | | | |
| 14 | 10/14/96 | 4.00 | 7.00 | 10.01 | 22.3 | 96.4 | CY |
| 15 | 10/15/96 | 3.99 | 7.00 | 10.01 | 22.2 | 98.8 | JEB |
| 16 | | | | | | | |
| 17 | 10/17/96 | 4.00 | 7.00 | 9.99 | 21.9 | 99.8 | CY |
| 18 | 10/18/96 | 4.00 | 7.00 | 10.01 | 22.4 | 97.8 | CY |
| 19 | | 4.00 | 7.00 | 10.01 | 22.4 | 97.8 | CY |
| 20 | | | | | | | |
| 21 | 10/21/96 | 4.00 | 7.00 | 10.05 | 22.4 | 92.912 | CY |
| 22 | 10/22/96 | 4.00 | 7.00 | 10.03 | 22.6 | 97.8 | CY |
| 23 | 10/23/96 | 4.00 | 7.00 | 9.95 | 21.9 | 98.6 | BSG |
| 24 | | | | | | | |
| 25 | | | | | | | |
| 26 | | | | | | | |
| 27 | | | | | | | |
| 28 | | | | | | | |
| 29 | | | | | | | |
| 30 | | | | | | | |
| 31 | | | | | | | |

Acceptance Limits =

 ± 0.05 pH UNITS ± 0.05 pH UNITS ± 0.05 pH UNITS

Refer to instructional manual if not within these limits