



TechLaw, Inc.
Environmental Services Assistance Team
16194 W. 45th Drive, Golden, CO 80403
303-312-7720

Task Order: 32
TDF: DG-220
LIMS: C101104
DCN #: EP8-5-5600
Contract: EP-W-06-33

Rico-Argentine – Surface Water – Nov 2010

Sample Identification
Cross-Reference Information

EPA Sample ID	ESAT LIMS ID
SLDRBG	C101104-01
SLDRBG	C101104-02
SLDRBG	C101104-03
SLDRBG DUP	C101104-04
SLDRBG DUP	C101104-05
SLDRBG DUP	C101104-06
SLDRMZ1a	C101104-07
SLDRMZ1a	C101104-08
SLDRMZ1a	C101104-09
SLDRMZ1b	C101104-10
SLDRMZ1b	C101104-11
SLDRMZ1b	C101104-12
SLDRMZ1c	C101104-13
SLDRMZ1c	C101104-14
SLDRMZ1c	C101104-15
SLDRMZ2	C101104-16
SLDRMZ2	C101104-17
SLDRMZ2	C101104-18
SLPO01	C101104-19
SLPO01	C101104-20
SLPO01	C101104-21
SLPO02	C101104-22
SLPO02	C101104-23
SLPO02	C101104-24
SLPO03	C101104-25
SLPO03	C101104-26
SLPO03	C101104-27

Rico-Argentina – Surface Water – Nov 2010

Page 2 of 3

Sample Identification Cross-Reference Information

EPA Sample ID	ESAT LIMS ID
SLPO04	C101104-28
SLPO04	C101104-29
SLPO04	C101104-30
SLPO05	C101104-31
SLPO05	C101104-32
SLPO05	C101104-33
SLSW01	C101104-34
SLSW01	C101104-35
SLSW01	C101104-36
SLSW02	C101104-37
SLSW02	C101104-38
SLSW02	C101104-39
SLSW03	C101104-40
SLSW03	C101104-41
SLSW03	C101104-42
SLSW04	C101104-43
SLSW04	C101104-44
SLSW04	C101104-45
SLSW05	C101104-46
SLSW05	C101104-47
SLSW05	C101104-48
SLSWDR3	C101104-49
SLSWDR3	C101104-50
SLSWDR3	C101104-51
SLSWDR4	C101104-52
SLSWDR4	C101104-53
SLSWDR4	C101104-54
SLSWDR6	C101104-55
SLSWDR6	C101104-56
SLSWDR6	C101104-57
SLSWDR7b	C101104-58
SLSWDR7b	C101104-59
SLSWDR7b	C101104-60
SLSWDR7c	C101104-61
SLSWDR7c	C101104-62
SLSWDR7c	C101104-63
SLSWFB	C101104-64
SLSWFB	C101104-65
SLSWFB	C101104-66
SLSWP06	C101104-67
SLSWP06	C101104-68
SLSWP06	C101104-69
SLSWP07a	C101104-70

Sample Identification
Cross-Reference Information

EPA Sample ID	ESAT LIMS ID
SLSWP07a	C101104-71
SLSWP07a	C101104-72
SLSWP07b	C101104-73
SLSWP07b	C101104-74
SLSWP07b	C101104-75
SLSWP08	C101104-76
SLSWP08	C101104-77
SLSWP08	C101104-78
SLSWP09	C101104-79
SLSWP09	C101104-80
SLSWP09	C101104-81
SLSWP10	C101104-82
SLSWP10	C101104-83
SLSWP10	C101104-84
SLSWP11	C101104-85
SLSWP11	C101104-86
SLSWP11	C101104-87
SLSWP12	C101104-88
SLSWP12	C101104-89
SLSWP12	C101104-90
SLSWP14	C101104-91
SLSWP14	C101104-92
SLSWP14	C101104-93
SLSWP15	C101104-94
SLSWP15	C101104-95
SLSWP15	C101104-96
SLSWPP	C101104-97
SLSWPP	C101104-98
SLSWPP	C101104-99

No: 11-17-10 Christner

Lab: ESAT - Region 8 Laboratory

Lab Phone: 303.312.7708

Contact Phone: 720.810.0795

Special Instructions:						SAMPLES TRANSFERRED FROM						F. D. [Signature]			
						CHAIN OF CUSTODY #						5CSE15D-sec			
Items/Reason		Relinquished by		Date	Received by	Date	Time	Items/Reason		Relinquished By		Date	Received by	Date	Time
3 coolers		M. J. [Signature]		11/19	[Signature]	11/19/10	14:10								

✓ 5CSE15D-sec

✓ 5CSE15D-sec

✓ 5CSE15D-sec

DateShipped: 11/18/2010
CarrierName: Hand Delivery
TDF No: DG-220

CHAIN OF CUSTODY RECORD

Site #: 11172010

Contact Name: Jan Christner

Contact Phone: 720.810.0795

No: 11-17-10 Christner

Lab: ESAT - Region 8 Laboratory



Lab Phone: 303.312.7708

Lab #	Sample #	Location	Analyses	Matrix	Collection Method	Collected	Sample Time	Numb Cont	Container	Preservative	MS/MSD
	SLPO01-111710	SLPO01	Total Metals	Surface Water	Grab	11/17/2010	11:19	1	1 L poly	none	
	SLPO01-111710	SLPO01	Dissolved Metals	Surface Water	Grab	11/17/2010	11:19	1	1 L poly	none	
	SLPO02-111710	SLPO02	Alkalinity	Surface Water	Grab	11/17/2010	10:48	1	250 mL poly	4 C	
	SLPO02-111710	SLPO02	Total Metals	Surface Water	Grab	11/17/2010	10:48	1	1 L poly	none	
	SLPO02-111710	SLPO02	Dissolved Metals	Surface Water	Grab	11/17/2010	10:48	1	1 L poly	none	
	SLPO03-111710	SLPO03	Alkalinity	Surface Water	Grab	11/17/2010	10:32	1	250 mL poly	4 C	
	SLPO03-111710	SLPO03	Total Metals	Surface Water	Grab	11/17/2010	10:32	1	1 L poly	none	
	SLPO03-111710	SLPO03	Dissolved Metals	Surface Water	Grab	11/17/2010	10:32	1	1 L poly	none	
	SLPO04-111710	SLPO04	Alkalinity	Surface Water	Grab	11/17/2010	10:15	1	250 mL poly	4 C	
	SLPO04-111710	SLPO04	Total Metals	Surface Water	Grab	11/17/2010	10:15	1	1 L poly	none	
	SLPO04-111710	SLPO04	Dissolved Metals	Surface Water	Grab	11/17/2010	10:15	1	1 L poly	none	
	SLPO05-111710	SLPO05	Alkalinity	Surface Water	Grab	11/17/2010	09:16	1	250 mL poly	4 C	
	SLPO05-111710	SLPO05	Total Metals	Surface Water	Grab	11/17/2010	09:16	1	1 L poly	none	
	SLPO05-111710	SLPO05	Dissolved Metals	Surface Water	Grab	11/17/2010	09:16	1	1 L poly	none	
	SLSE02-111710	SLSE02	Metals	Sediment	Grab	11/17/2010	09:35	1	8 oz glass	4 C	W/leaker
	SLSE04-111710	SLSE04	Metals	Sediment	Grab	11/17/2010	09:41	1	8 oz glass	4 C	10/10/05
	SLSE05-111710	SLSE05	Metals	Sediment	Grab	11/17/2010	09:45	1	8 oz glass	4 C	
	SLSE10-111610	SLSE10	Metals	Sediment	Grab	11/16/2010	15:30	1	8 oz glass	4 C	
	SLSE15-06-111610	SLSE15-06	Metals	Sediment	Grab	11/16/2010	15:20	1	8 oz glass	4 C	

Special Instructions:

SAMPLES TRANSFERRED FROM	
1	2
3	4
5	6
7	8
9	10
11	12
13	14
15	16
17	18
19	20
21	22
23	24
25	26
27	28
29	30
31	32
33	34
35	36
37	38
39	40
41	42
43	44
45	46
47	48
49	50
51	52
53	54
55	56
57	58
59	60
61	62
63	64
65	66
67	68
69	70
71	72
73	74
75	76
77	78
79	80
81	82
83	84
85	86
87	88
89	90
91	92
93	94
95	96
97	98
99	100

CHAIN OF CUSTODY #

Items/Reason	Relinquished by	Date	Received by	Date	Time	Items/Reason	Relinquished By	Date	Received by	Date	Time
				11/19/10	14:10						

• Site #: 11172010

Contact Phone: 720.810.0795

Lab Phone: 303.312.7708

Items/Reason	Relinquished by	Date	Received by	Date	Time
	[Signature]		[Signature]	11/18/10	14:10

Site #: 11172010

Contact Phone: 720.810.0795

Lab Phone: 303.312.7708

Lab #	Sample #	Location	Analyses	Matrix	Collection Method	Collected	Sample Time	Number Cont	Container	Preservative	MS/MSD
	SLSW05-111710	SLSW05	Alkalinity	Surface Water	Grab	11/17/2010	11:55	1	250 mL poly	4 C	
	SLSW05-111710	SLSW05	Total Metals	Surface Water	Grab	11/17/2010	11:55	1	1 L poly	none	
	SLSW05-111710	SLSW05	Dissolved Metals	Surface Water	Grab	11/17/2010	11:55	1	1 L poly	none	
	SLSWDR3-111610	SLSWDR3	Alkalinity	Surface Water	Grab	11/16/2010	09:20	1	250 mL poly	4 C	
	SLSWDR3-111610	SLSWDR3	Total Metals	Surface Water	Grab	11/16/2010	09:20	1	1 L poly	HNO3 pH<2	
	SLSWDR3-111610	SLSWDR3	Dissolved Metals	Surface Water	Grab	11/16/2010	09:20	1	1 L poly	HNO3 pH<2	
	SLSWDR4-111610	SLSWDR4	Alkalinity	Surface Water	Grab	11/16/2010	09:30	1	250 mL poly	4 C	
	SLSWDR4-111610	SLSWDR4	Total Metals	Surface Water	Grab	11/16/2010	09:30	1	1 L poly	HNO3 pH<2	
	SLSWDR4-111610	SLSWDR4	Dissolved Metals	Surface Water	Grab	11/16/2010	09:30	1	1 L poly	HNO3 pH<2	
	SLSWDR6-111610	SLSWDR6	Alkalinity	Surface Water	Grab	11/16/2010	14:00	1	250 mL poly	4 C	
	SLSWDR6-111610	SLSWDR6	Total Metals	Surface Water	Grab	11/16/2010	14:00	1	1 L poly	HNO3 pH<2	
	SLSWDR6-111610	SLSWDR6	Dissolved Metals	Surface Water	Grab	11/16/2010	14:00	1	1 L poly	HNO3 pH<2	
	SLSWDR7b-111610	SLSWDR7b	Alkalinity	Surface Water	Grab	11/16/2010	16:50	1	250 mL poly	4 C	
	SLSWDR7b-111610	SLSWDR7b	Total Metals	Surface Water	Grab	11/16/2010	16:50	1	1 L poly	HNO3 pH<2	
	SLSWDR7b-111610	SLSWDR7b	Dissolved Metals	Surface Water	Grab	11/16/2010	16:50	1	1 L poly	HNO3 pH<2	
	SLSWDR7c-111610	SLSWDR7c	Alkalinity	Surface Water	Grab	11/16/2010	16:50	1	250 mL poly	4 C	
	SLSWDR7c-111610	SLSWDR7c	Total Metals	Surface Water	Grab	11/16/2010	16:50	1	1 L poly	HNO3 pH<2	
	SLSWDR7c-111610	SLSWDR7c	Dissolved Metals	Surface Water	Grab	11/16/2010	16:50	1	1 L poly	HNO3 pH<2	
	SLSWFB-111710	SLSWFB	Alkalinity	Surface Water	Grab	11/17/2010	18:00	1	250 mL poly	4 C	

Special Instructions:	SAMPLES TRANSFERRED FROM
	CHAIN OF CUSTODY #

[illegible]

DateShipped: 11/18/2010
CarrierName: Hand Delivery
TDF No: DG-220

CHAIN OF CUSTODY RECORD

Site #: 11172010
Contact Name: Jan Christner
Contact Phone: 720.810.0795

No: 11-17-10 Christner

Lab: ESAT - Region 8 Laboratory
Lab Phone: 303.312.7708

Lab #	Sample #	Location	Analyses	Matrix	Collection Method	Collected	Sample Time	Number of Containers	Container	Preservative	MS/MSD
	SLSWFB-111710	SLSWFB	Total Metals	Surface Water	Grab	11/17/2010	18:00	1	250 mL poly	none	
	SLSWFB-111710	SLSWFB	Dissolved Metals	Surface Water	Grab	11/17/2010	18:00	1	250 mL poly	none	
	SLSWP06-111610	SLSWP06	Alkalinity	Surface Water	Grab	11/16/2010	13:45	1	250 mL poly	4 C	
	SLSWP06-111610	SLSWP06	Total Metals	Surface Water	Grab	11/16/2010	13:45	1	1 L poly	HNO3 pH<2	
	SLSWP06-111610	SLSWP06	Dissolved Metals	Surface Water	Grab	11/16/2010	13:45	1	1 L poly	HNO3 pH<2	
	SLSWP07a-111610	SLSWP07a	Alkalinity	Surface Water	Grab	11/16/2010	12:15	1	250 mL poly	4 C	
	SLSWP07a-111610	SLSWP07a	Total Metals	Surface Water	Grab	11/16/2010	12:15	1	1 L poly	HNO3 pH<2	
	SLSWP07a-111610	SLSWP07a	Dissolved Metals	Surface Water	Grab	11/16/2010	12:15	1	1 L poly	HNO3 pH<2	
	SLSWP07b-111610	SLSWP07b	Alkalinity	Surface Water	Grab	11/16/2010	12:15	1	250 mL poly	4 C	
	SLSWP07b-111610	SLSWP07b	Total Metals	Surface Water	Grab	11/16/2010	12:15	1	1 L poly	HNO3 pH<2	
	SLSWP07b-111610	SLSWP07b	Dissolved Metals	Surface Water	Grab	11/16/2010	12:15	1	1 L poly	HNO3 pH<2	
	SLSWP08-111610	SLSWP08	Alkalinity	Surface Water	Grab	11/16/2010	12:00	1	250 mL poly	4 C	
	SLSWP08-111610	SLSWP08	Total Metals	Surface Water	Grab	11/16/2010	12:00	1	1 L poly	HNO3 pH<2	
	SLSWP08-111610	SLSWP08	Dissolved Metals	Surface Water	Grab	11/16/2010	12:00	1	1 L poly	HNO3 pH<2	
	SLSWP09-111610	SLSWP09	Alkalinity	Surface Water	Grab	11/16/2010	11:50	1	250 mL poly	4 C	
	SLSWP09-111610	SLSWP09	Total Metals	Surface Water	Grab	11/16/2010	11:50	1	1 L poly	HNO3 pH<2	
	SLSWP09-111610	SLSWP09	Dissolved Metals	Surface Water	Grab	11/16/2010	11:50	1	1 L poly	HNO3 pH<2	
	SLSWP10-111610	SLSWP10	Alkalinity	Surface Water	Grab	11/16/2010	11:30	1	250 mL poly	4 C	
	SLSWP10-111610	SLSWP10	Total Metals	Surface Water	Grab	11/16/2010	11:30	1	1 L poly	HNO3 pH<2	

Special Instructions:	SAMPLES TRANSFERRED FROM
	CHAIN OF CUSTODY #

[illegible]

Contact Phone: 720.810.0795

Lab Phone: 303.312.7708

[illegible]

ESAT Technical Direction Form

Contract No. EPW06033
EPA Region 8

C101104

Site ID: 08BU
TDF ID: DG-220

Date Issued: 11/16/2010
Date Updated:

Date
Closed By:

Name: Rico - Argentine Samples

Details: The Contractor shall analyze several soil and aqueous samples collected by UOS from the Rico-Argentine Superfund Site. The samples are expected to be delivered to the ESAT R8 Laboratory on approximately 11/19/10. The samples are expected to be analyzed for metals and alkalinity as indicated on the COCs and in the analytical information section below.

The OSC for the site is Steve Way.

Analytical

MATRIX

☒ Water ☒ Soils ☐ Vegetation ☐ Biota

WET CHEM

☐ TSS ☐ TDS ☐ DOC ☒ Alk ☐ Chloride ☐ Sulfate ☐ Fluoride ☐ Nitrate ☐ Nitrite
Other

METALS

☒ Dissolved ☒ Total Recoverable ☐ Total ☒ Hardness (Calc) ✓
200.7: ☐ Ag ☒ Al ☐ As ☒ Ba ☐ Be ☐ B ☒ Ca ☐ Cd ☐ Co ☐ Cr ☐ Cu ☒ Fe ☒ K ☒ Mg
☐ Mn ☐ Mo ☒ Na ☐ Ni ☐ Pb ☐ Sb ☐ Se ☐ Sr ☐ Ti ☐ Tl ☐ V ☒ Zn ☐ SiO2
200.8: ☒ Ag ☐ Al ☒ As ☐ Ba ☐ Be ☒ Cd ☒ Co ☒ Cr ☒ Cu ☒ Mn ☐ Mo ☒ Ni ☒ Pb ☐ Sb
☒ Se ☐ Th ☐ Tl ☐ U ☐ V ☐ Zn
7470/7471/747 ☐ Hg

FIBERS

☐ PLM ☐ TEM

Deliverable

ID	Description	Due Date	Submission Date
1	Provide final deliverable package to Task Monitor no later than 30 days after delivery of samples.	12/20/2010	



TechLaw, Inc.
Environmental Services Assistance Team
16194 W. 45th Drive, Golden, CO 80403
303-312-7720

Task Order: 32
Contract: EP-W-06-33
Valid: April 2010 – May 2011

Certificates of Analysis

Valid through May 1, 2010

Perkin Elmer Optima ICP-OE

Perkin Elmer ELAN 6000 ICP-MS

NIPPON NIC MA2000

Perkin Elmer FIMS 100

- Initial Calibration Verification (ICV) Standards
- Laboratory Check Standards (LCS)
- Matrix Spike Solutions
- Interference Check (ICSA / AB) Standards

13594
4/6/10

tel: 800.669.6799 · 540.585.3030

fax: 540.585.3012

info@inorganicventures.com

- 1.0 INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



2.0 DESCRIPTION OF CRM 1000 µg/mL Aluminum in 3% (v/v) HNO₃

Catalog Number: CGAL1-1, CGAL1-2, and CGAL1-5

Lot Number: C2-AL04078

Starting Material: Al ingot

Starting Material Purity (%): 99.998788

Starting Material Lot No: C14S012

Matrix: 3% (v/v) HNO₃

EXPIRES
1-May-2011

3.0 CERTIFIED VALUES AND UNCERTAINTIES

Certified Concentration: 1003 ± 5 µg/mL

Certified Density: 1.017 g/mL (measured at 22° C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

"Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

4.1 Assay Method #1 1003 ± 5 µg/mL

ICP Assay NIST SRM 3101a Lot Number: 060502

Assay Method #2 1002 ± 8 µg/mL

EDTA NIST SRM 928 Lot Number: 928

13596
4/6/10

tel: 800.669.6799 · 540.585.3030

fax: 540.585.3012

info@inorganicventures.com

- 1.0 **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0 **DESCRIPTION OF CRM** **1000 µg/mL Antimony in 1% (v/v) HNO₃ / 3% Tartaric Acid**
- Catalog Number: CGSB1-1, CGSB1-2 and CGSB1-5
- Lot Number: **C2-SB02121**
- Starting Material: Sb shot
- Starting Material Purity (%): 99.996681
- Starting Material Lot No: R1105SBA1
- Matrix: 1% (v/v) HNO₃ / 3% Tartaric Acid

EXPIRES
1-May-2011

3.0 **CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 990 ± 4 µg/mL

Certified Density: 1.021 g/mL (measured at 22° C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of $k = 2$.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i^2)^{1/2}]}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 **TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS**

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

4.1 **Assay Method #1** **990 ± 4 µg/mL (avg. of 2 runs)**

ICP Assay NIST SRM 3102A Lot Number: 061229

Assay Method #2 **1000 ± 5 µg/mL**

Calculated NIST SRM Lot Number: See Sec. 4.2

13597
4/6/10

tel: 800.669.6799 · 540.585.3030

fax: 540.585.3012

info@inorganicventures.com

- 1.0 INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0 DESCRIPTION OF CRM** **1000 µg/mL Arsenic in 2% (v/v) HNO₃**
- Catalog Number: CGAS1-1, CGAS1-2, and CGAS1-5
- Lot Number: **C2-AS02061**
- Starting Material: As pieces
- Starting Material Purity (%): 99.999032
- Starting Material Lot No: R1107ASB1
- Matrix: 2% (v/v) HNO₃

EXPIRES**1-May-2011**

3.0 CERTIFIED VALUES AND UNCERTAINTIES

Certified Concentration: 999 ± 4 µg/mL

Certified Density: 1.011 g/mL (measured at 20 ± 1°C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of $k = 2$.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- 4.1 Assay Method #1** **999 ± 4 µg/mL**
ICP Assay NIST SRM 3103a Lot Number: 010713
- Assay Method #2** **1001 ± 5 µg/mL**
Calculated NIST SRM Lot Number: See Sec. 4.2

- 1.0 **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0 **DESCRIPTION OF CRM** **1000 µg/mL Barium in tr. HNO₃**
- Catalog Number: CGBA1-1, CGBA1-2, and CGBA1-5
- Lot Number: **C2-BA02050**
- Starting Material: Ba(NO₃)₂
- Starting Material Purity (%): 99.999911
- Starting Material Lot No: W500A
- Matrix: tr. HNO₃

EXPIRES
1-May-2011

3.0 **CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 1,000 ± 2 µg/mL

Certified Density: 1.000 g/mL (measured at 22° C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i^2)^{1/2}]}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 **TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS**

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- 4.1 **Assay Method #1** **999 ± 4 µg/mL**
ICP Assay NIST SRM 3104a Lot Number: 070222
- Assay Method #2** **1,000 ± 2 µg/mL**
Gravimetric NIST SRM Lot Number: See Sec. 4.2

13600 13589
4/6/10

- 1.0 **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0 **DESCRIPTION OF CRM** **1000 µg/mL Beryllium in 3% (v/v) HNO₃**

Catalog Number: CGBE1-1, CGBE1-2, and CGBE1-5

Lot Number: **C2-BE01124**

Starting Material: Be(OOCCH₃)₂

Starting Material Purity (%): 99.999948

Starting Material Lot No: 0801-1

Matrix: 3% (v/v) HNO₃

EXPIRES
1-May-2011

- 3.0 **CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 1003 ± 4 µg/mL

Certified Density: 1.022 g/mL (measured at 22° C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

- 4.0 **TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS**

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- 4.1 **Assay Method #1** **1003 ± 4 µg/mL**

ICP Assay NIST SRM 3105a Lot Number: 892707

- Assay Method #2** **1001 ± 5 µg/mL**

Calculated NIST SRM Lot Number: See Sec. 4.2

13589
4/6/10

- 1.0** **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0** **DESCRIPTION OF CRM** **1000 µg/mL Boron in H₂O**
- Catalog Number: CGB1-1, CGB1-2, and CGB1-5
- Lot Number: **C2-B02083**
- Starting Material: H₃BO₃
- Starting Material Purity (%): 99.999998
- Starting Material Lot No: OV0133
- Matrix: H₂O

EXPIRES
1-May-2011

3.0 **CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 998 ± 3 µg/mL

Certified Density: 1.000 g/mL (measured at 22° C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of $k = 2$.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 **TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS**

• "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

• This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- 4.1** **Assay Method #1** **998 ± 3 µg/mL**
ICP Assay NIST SRM 3107 Lot Number: 070514
- Assay Method #2** **1,001 ± 5 µg/mL**
Calculated NIST SRM Lot Number: See Sec. 4.2

13601
4/6/10tel: 800.669.6799 · 540.585.3030
fax: 540.585.3012
info@inorganicventures.com

- 1.0 **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



2.0 **DESCRIPTION OF CRM** **10000 µg/mL Bromide (NH₄) in H₂O**

Catalog Number: CGICBR10-1 and CGICBR10-5

Lot Number: **A2-BR01066**

Starting Material: NH₄Br

Starting Material Purity (%): 99.998491

Starting Material Lot No: DI05205EU

Matrix: H₂O

EXPIRES
1-May-2011

3.0 **CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 9993 ± 34 µg/mL

Certified Density: 1.004 g/mL (measured at 20 ± 1°C)

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume and the fixed error reported on the NIST SRM certificate of analysis)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of $k = 2$.

4.0 **TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS**

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- 4.1 **Assay Method #1** **9993 ± 34 µg/mL**
IC Assay NIST SRM 3184 Lot Number: 020701
- Assay Method #2** **9976 ± 47 µg/mL**
Volhard NIST SRM 999b Lot Number: 999b

13502
4/6/10

- 1.0 INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0 DESCRIPTION OF CRM** **1000 µg/mL Cadmium in 3% (v/v) HNO₃**
- Catalog Number: CGCD1-1, CGCD1-2, and CGCD1-5
- Lot Number: **C2-CD02021**
- Starting Material: Cd shot
- Starting Material Purity (%): 99.999656
- Starting Material Lot No: R1205CDA1
- Matrix: 3% (v/v) HNO₃

EXPIRES
1-May-2011

3.0 CERTIFIED VALUES AND UNCERTAINTIES

Certified Concentration: 998 ± 3 µg/mL

Certified Density: 1.015 g/mL (measured at 20 ± 1°C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- 4.1 Assay Method #1 995 ± 3 µg/mL**
ICP Assay NIST SRM 3108 Lot Number: 060531
- Assay Method #2 998 ± 3 µg/mL**
EDTA NIST SRM 928 Lot Number: 928

13603
466/10

- 1.0 INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0 DESCRIPTION OF CRM** **1000 µg/mL Calcium in 0.1% (v/v) HNO₃**
- Catalog Number: CGCA1-1, CGCA1-2, and CGCA1-5
- Lot Number: **C2-CA03123**
- Starting Material: CaCO₃
- Starting Material Purity (%): 99.998463
- Starting Material Lot No: C808CAA1
- Matrix: 0.1% (v/v) HNO₃

EXPIRES
1-May-2011

3.0 CERTIFIED VALUES AND UNCERTAINTIES

Certified Concentration: 997 ± 3 µg/mL

Certified Density: 1.001 g/mL (measured at 20 ± 1°C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of $k = 2$.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- 4.1 Assay Method #1 993 ± 3 µg/mL**
ICP Assay NIST SRM 3109a Lot Number: 050825
- Assay Method #2 997 ± 3 µg/mL**
EDTA NIST SRM 928 Lot Number: 928

132605
4/6/10

- 1.0** **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0** **DESCRIPTION OF CRM** **1000 µg/mL Chromium (+3) in 2% (v/v) HNO₃**
- Catalog Number: CGCR(3)1-1, CGCR(3)1-2, and CGCR(3)1-5
- Lot Number: **C2-CR03027**
- Starting Material: Cr pieces
- Starting Material Purity (%): 99.993508
- Starting Material Lot No: R800A
- Matrix: 2% (v/v) HNO₃

EXPIRES
1-May-2011

3.0 **CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 990 ± 3 µg/mL

Certified Density: 1.013 g/mL (measured at 22° C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of $k = 2$.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 **TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS**

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- 4.1** **Assay Method #1** **990 ± 3 µg/mL**
ICP Assay NIST SRM 3112a Lot Number: 030730
- Assay Method #2** **1000 ± 5 µg/mL**
Calculated NIST SRM Lot Number: See Sec. 4.2

13686
4/6/10

- 1.0** **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."

2.0 DESCRIPTION OF CRM **1000 µg/mL Cobalt in 3% (v/v) HNO₃**

Catalog Number: CGCO1-1, CGCO1-2, and CGCO1-5

Lot Number: **C2-CO02022**

Starting Material: Co powder

Starting Material Purity (%): 99.997920

Starting Material Lot No: PW407COA1

Matrix: 3% (v/v) HNO₃

EXPIRES
1-May-2011

3.0 CERTIFIED VALUES AND UNCERTAINTIES

Certified Concentration: 995 ± 3 µg/mL

Certified Density: 1.018 g/mL (measured at 20 ± 1°C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of $k = 2$.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term "in-house std." is specified.

4.1 Assay Method #1 **995 ± 3 µg/mL**

ICP Assay NIST SRM 3113 Lot Number: 00630

Assay Method #2 **1,004 ± 4 µg/mL**

EDTA NIST SRM 928 Lot Number: 928

13807
4/6/10

- 1.0** **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0** **DESCRIPTION OF CRM** **1000 µg/mL Copper in 3% (v/v) HNO₃**
- Catalog Number: CGCU1-1, CGCU1-2, and CGCU1-5
- Lot Number: **C2-CU02116**
- Starting Material: Cu shot
- Starting Material Purity (%): 99.999779
- Starting Material Lot No: R508CUA2
- Matrix: 3% (v/v) HNO₃

EXPIRES
1-May-2011

3.0 **CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 992 ± 3 µg/mL

Certified Density: 1.015 g/mL (measured at 20 ± 1°C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 **TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS**

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- 4.1** **Assay Method #1** **994 ± 4 µg/mL**
ICP Assay NIST SRM 3114 Lot Number: 011017
- Assay Method #2** **992 ± 3 µg/mL**
EDTA NIST SRM 928 Lot Number: 928

13608
466/100

- 1.0** **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0** **DESCRIPTION OF CRM** **1000 µg/mL Gadolinium in 7% (v/v) HNO₃**
- Catalog Number: CGGD1-1, CGGD1-2, and CGGD1-5
- Lot Number: **C2-GD01038**
- Starting Material: Gd₂O₃
- Starting Material Purity (%): 99.999400
- Starting Material Lot No: GD-0-5-026
- Matrix: 7% (v/v) HNO₃

EXPIRES
1-May-2011

3.0 **CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 1,001 ± 3 µg/mL

Certified Density: 1.035 g/mL (measured at 22° C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 **TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS**

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- 4.1** **Assay Method #1** **1,001 ± 3 µg/mL**
 EDTA NIST SRM 928 Lot Number: 928
- Assay Method #2** **999 ± 3 µg/mL**
 ICP Assay NIST SRM 3118a Lot Number: 992004

13609
4/6/10

- 1.0 **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0 **DESCRIPTION OF CRM** **1000 µg/mL Iron in 2% (v/v) HNO₃**
Catalog Number: CGFE1-1, CGFE1-2, and CGFE1-5
Lot Number: **C2-FE03104**
Starting Material: Fe powder
Starting Material Purity (%): 99.997283
Starting Material Lot No: R1207FEA1
Matrix: 2% (v/v) HNO₃

EXPIRES
1-May-2011

3.0 **CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 996 ± 3 µg/mL

Certified Density: 1.011 g/mL (measured at 20 ± 1°C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of $k = 2$.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 **TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS**

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- 4.1 **Assay Method #1** **993 ± 4 µg/mL**
ICP Assay NIST SRM 3126a Lot Number: 051031
- Assay Method #2** **996 ± 3 µg/mL**
EDTA NIST SRM 928 Lot Number: 928

- 1.0 INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0 DESCRIPTION OF CRM** **1000 µg/mL Lead in 0.5% (v/v) HNO₃**
- Catalog Number: CGPB1-1, CGPB1-2, and CGPB1-5
- Lot Number: **C2-PB03013**
- Starting Material: Pb(NO₃)₂
- Starting Material Purity (%): 99.999554
- Starting Material Lot No: E1007PBA1
- Matrix: 0.5% (v/v) HNO₃

EXPIRES
1-May-2011

3.0 CERTIFIED VALUES AND UNCERTAINTIES

Certified Concentration: 1004 ± 4 µg/mL

Certified Density: 1.001 g/mL (measured at 20 ± 1°C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors (Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- 4.1 Assay Method #1 999 ± 5 µg/mL**
ICP Assay NIST SRM 3128 Lot Number: 030721
- Assay Method #2 1004 ± 4 µg/mL**
EDTA NIST SRM 928 Lot Number: 928

- 1.0 **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0 **DESCRIPTION OF CRM** **1000 µg/mL Lithium in 0.1% (v/v) HNO₃**

Catalog Number: CGLI1-1, CGLI1-2, and CGLI1-5

Lot Number: **C2-LI02113**

Starting Material: Li₂CO₃

Starting Material Purity (%): 99.997165

Starting Material Lot No: 1123

Matrix: 0.1% (v/v) HNO₃

EXPIRES
1-May-2011

- 3.0 **CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 997 ± 1 µg/mL

Certified Density: 1.005 g/mL (measured at 22° C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

- 4.0 **TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS**

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- 4.1 **Assay Method #1** **994 ± 3 µg/mL**

ICP Assay NIST SRM 3129a Lot Number: 000505

- Assay Method #2** **997 ± 1 µg/mL**

Gravimetric NIST SRM Lot Number: See Sec. 4.2

136134
4/6/10

- 1.0 INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0 DESCRIPTION OF CRM** **1000 µg/mL Magnesium in 0.1% (v/v) HNO₃**
- Catalog Number: CGMG1-1, CGMG1-2, and CGMG1-5
- Lot Number: **C2-MG03084**
- Starting Material: Mg metal
- Starting Material Purity (%): 99.999531
- Starting Material Lot No: 00265067003
- Matrix: 0.1% (v/v) HNO₃

EXPIRES
1-May-2011

3.0 CERTIFIED VALUES AND UNCERTAINTIES

Certified Concentration: 997 ± 3 µg/mL

Certified Density: 1.002 g/mL (measured at 20 ± 1°C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

"Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- 4.1 Assay Method #1 995 ± 3 µg/mL**
ICP Assay NIST SRM 3131a Lot Number: 050302
- Assay Method #2 997 ± 3 µg/mL**
EDTA NIST SRM 928 Lot Number: 928

13615
4/6/10

- 1.0 INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0 DESCRIPTION OF CRM** **1000 µg/mL Manganese in 3% HNO₃ (v/v)**
- Catalog Number: CGMN1-1, CGMN1-2, and CGMN1-5
- Lot Number: **B2-MN02070**
- Starting Material: Mn pieces
- Starting Material Purity (%): 99.994492
- Starting Material Lot No: R806MNA1
- Matrix: 3% HNO₃ (v/v)

EXPIRES
1-May-2011**3.0 CERTIFIED VALUES AND UNCERTAINTIES****Certified Concentration:** 1,005 ± 2 µg/mL**Certified Density:** 1.016 g/mL (measured at 22° C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of $k = 2$.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

 (\bar{x}) = mean x_i = individual results n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

 $\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

"Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- 4.1 Assay Method #1** **1,005 ± 2 µg/mL**
ICP Assay NIST SRM 3132 Lot Number: 050429
- Assay Method #2** **1,004 ± 2 µg/mL**
EDTA NIST SRM 928 Lot Number: 928

13616
4/6/10

- 1.0 **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0 **DESCRIPTION OF CRM** **1000 µg/mL Mercury in 5% (v/v) HNO₃**
- Catalog Number: CGHG1-1, CGHG1-2, and CGHG1-5
- Lot Number: **C2-HG02070**
- Starting Material: Hg metal
- Starting Material Purity (%): 99.999792
- Starting Material Lot No: R307HGA1
- Matrix: 5% (v/v) HNO₃

EXPIRES
1-May-2011

3.0 **CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 1003 ± 3 µg/mL

Certified Density: 1.027 g/mL (measured at 22° C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 **TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS**

"Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- 4.1 **Assay Method #1** **999 ± 5 µg/mL**
ICP Assay NIST SRM 3133 Lot Number: 061204

- Assay Method #2** **1003 ± 3 µg/mL**
EDTA NIST SRM 928 Lot Number: 928

13617
4/6/10

- 1.0** **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0** **DESCRIPTION OF CRM** **1000 µg/mL Molybdenum in H₂O / tr. NH₄OH**
- Catalog Number: CGMO1-1, CGMO1-2, and CGMO1-5
- Lot Number: **C2-MO02032**
- Starting Material: (NH₄)₆Mo₇O₂₄·xH₂O
- Starting Material Purity (%): 99.998755
- Starting Material Lot No: P704MOA1
- Matrix: H₂O / tr. NH₄OH

EXPIRES
1-May-2011

3.0 **CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 999 ± 3 µg/mL

Certified Density: 0.999 g/mL (measured at 22° C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 **TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS**

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- 4.1** **Assay Method #1** **999 ± 3 µg/mL (Avg. of 2 Runs)**
ICP Assay NIST SRM 3134 Lot Number: 891307
- Assay Method #2** **1001 ± 5 µg/mL**
Calculated NIST SRM Lot Number: See Sec. 4.2

- 1.0 **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0 **DESCRIPTION OF CRM** **1000 µg/mL Nickel in 2% (v/v) HNO₃**
- Catalog Number: CGNI1-1, CGNI1-2, and CGNI1-5
- Lot Number: **C2-NI02062**
- Starting Material: Ni pieces
- Starting Material Purity (%): 99.999033
- Starting Material Lot No: E25T014
- Matrix: 2% (v/v) HNO₃

13618
4/6/10**EXPIRES**
1-May-2011

3.0 **CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 1003 ± 3 µg/mL

Certified Density: 1.012 g/mL (measured at 22° C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)


4.0 **TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS**

• "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

• This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- 4.1 **Assay Method #1** **998 ± 3 µg/mL**
ICP Assay NIST SRM 3136 Lot Number: 000612
- Assay Method #2** **1003 ± 3 µg/mL**
EDTA NIST SRM 928 Lot Number: 928

- 1.0 **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



2.0 **DESCRIPTION OF CRM** **1000 µg/mL Potassium in 1% HNO₃**

Catalog Number: CGK1-1, CGK1-2, and CGK1-5

Lot Number: **C2-K03004**

Starting Material: KNO₃

Starting Material Purity (%): 99.996911

Starting Material Lot No: B19P01

Matrix: 1% HNO₃

EXPIRES
1-May-2011

3.0 **CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 1006 ± 4 µg/mL

Certified Density: 1.004 g/mL (measured at 22° C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 **TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS**

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

4.1 **Assay Method #1** **1006 ± 4 µg/mL**

ICP Assay NIST SRM 3141a Lot Number: 051220

Assay Method #2 **1007 ± 2 µg/mL**

Gravimetric NIST SRM Lot Number: See Sec. 4.2

13624
4/6/10

- 1.0 INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025, "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0 DESCRIPTION OF CRM** **1000 µg/mL Selenium(+4) in 2% (v/v) HNO₃**
- Catalog Number: CGSE(4)1-1, CGSE(4)1-2, and CGSE(4)1-5
- Lot Number: **C2-SE02010**
- Starting Material: Se shot
- Starting Material Purity (%): 99.999239
- Starting Material Lot No: B14A38
- Matrix: 2% (v/v) HNO₃

EXPIRES
1-May-2011

3.0 CERTIFIED VALUES AND UNCERTAINTIES

Certified Concentration: 992 ± 4 µg/mL

Certified Density: 1.010 g/mL (measured at 20 ± 1°C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i^2)]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume and the fixed error reported on the NIST SRM certificate of analysis)

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- 4.1 Assay Method #1** **992 ± 4 µg/mL (avg. of 2 runs)**
ICP Assay NIST SRM 3149 Lot Number: 992106
- Assay Method #2** **1000 ± 5 µg/mL**
Calculated NIST SRM Lot Number: See Sec. 4.2

13622
H6/10

- 1.0 INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0 DESCRIPTION OF CRM** **1000 µg/mL Silica in 1% (v/v) HNO₃ / tr. HF**
- Catalog Number: CGSIO1-1, CGSIO1-2, and CGSIO1-5
- Lot Number: **C2-SI02130**
- Starting Material: SiO₂
- Starting Material Purity (%): 99.999846
- Starting Material Lot No: 1015122300.00
- Matrix: 1% (v/v) HNO₃ / tr. HF

EXPIRES
1-May-2011**3.0 CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 999 ± 2 µg/mL

Certified Density: 1.005 g/mL (measured at 20 ± 1°C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i^2)^{1/2}]}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

4.1 Assay Method #1 **999 ± 2 µg/mL (Avg. of 2 Runs)**

ICP Assay NIST SRM 3150 Lot Number: 071204

13023
4/6/10

- 1.0** **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0** **DESCRIPTION OF CRM** **1000 µg/mL Silver in 5% (v/v) HNO₃**
- Catalog Number: CGAG1-1, CGAG1-2, and CGAG1-5
- Lot Number: **C2-AG02042**
- Starting Material: Ag shot
- Starting Material Purity (%): 99.991024
- Starting Material Lot No: E308AGA1
- Matrix: 5% (v/v) HNO₃

EXPIRES
1-May-2011

3.0 **CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 1001 ± 1 µg/mL

Certified Density: 1.024 g/mL (measured at 20 ± 1°C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i^2)^{1/2}]}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume and the fixed error reported on the NIST SRM certificate of analysis)

4.0 **TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS**

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- 4.1** **Assay Method #1** **996 ± 5 µg/mL**
ICP Assay NIST SRM 3151 Lot Number: 992212
- Assay Method #2** **1001 ± 1 µg/mL**
Volhard NIST SRM 999b Lot Number: 999b

17624
4/6/10

- 1.0 **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0 **DESCRIPTION OF CRM** **1000 µg/mL Sodium in 0.1% (v/v) HNO₃**

Catalog Number: CGNA1-1, CGNA1-2, and CGNA1-5

Lot Number: **C2-NA03074**

Starting Material: Na₂CO₃

Starting Material Purity (%): 99.997122

Starting Material Lot No: C18157

Matrix: 0.1% (v/v) HNO₃

EXPIRES
1-May-2011

- 3.0 **CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 1,006 ± 5 µg/mL

Certified Density: 1.001 g/mL (measured at 20 ± 1°C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of $k = 2$.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i^2)^{1/2}]}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

- 4.0 **TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS**

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- 4.1 **Assay Method #1** **1,006 ± 5 µg/mL**
ICP Assay NIST SRM 3152a Lot Number: 010728

- Assay Method #2** **1,008 ± 5 µg/mL**
Gravimetric NIST SRM Lot Number: See Sec. 4.2

- 1.0 **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0 **DESCRIPTION OF CRM** **1000 µg/mL Strontium in 0.1% (v/v) HNO₃**

Catalog Number: CGSR1-1, CGSR1-2, and CGSR1-5

Lot Number: **C2-SR02024**

Starting Material: SrCO₃

Starting Material Purity (%): 99.998364

Starting Material Lot No: W999A

Matrix: 0.1% (v/v) HNO₃

EXPIRES
1-May-2011

- 3.0 **CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 1002 ± 3 µg/mL

Certified Density: 1.000 g/mL (measured at 22° C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

- 4.0 **TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS**

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- 4.1 **Assay Method #1** **1002 ± 3 µg/mL**

EDTA NIST SRM 928 Lot Number: 928

- Assay Method #2** **1000 ± 5 µg/mL**

ICP Assay NIST SRM 3153a Lot Number: 990906

- 1.0** **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0** **DESCRIPTION OF CRM** **1000 µg/mL Thallium in 0.7% (v/v) HNO₃**
- Catalog Number: CGTL1-1, CGTL1-2, and CGTL1-5
- Lot Number: **D2-TL01111**
- Starting Material: TINO₃
- Starting Material Purity (%): 99.999671
- Starting Material Lot No: 1479
- Matrix: 0.7% (v/v) HNO₃

EXPIRES
1-May-2011

3.0 **CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 1,000 ± 4 µg/mL

Certified Density: 1.003 g/mL (measured at 20 ± 1°C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of $k = 2$.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i^2)]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 **TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS**

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- 4.1** **Assay Method #1** **1,000 ± 4 µg/mL (avg. of 2 runs)**
ICP Assay NIST SRM 3158 Lot Number: 993012
- Assay Method #2** **1,001 ± 5 µg/mL**
Calculated NIST SRM Lot Number: See Sec. 4.2

1.0 **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



2.0 **DESCRIPTION OF CRM** **1000 µg/mL Thorium in 4% (v/v) HNO₃**

Catalog Number: CGTH1-1, CGTH1-2, and CGTH1-5

Lot Number: **B2-TH01073**

Starting Material: Th(NO₃)₄•4H₂O

Starting Material Purity (%): 99.998809

Starting Material Lot No: X-25828-7

Matrix: 4% (v/v) HNO₃

EXPIRES

1-May-2011

3.0 **CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 1,003 ± 3 µg/mL

Certified Density: 1.021 g/mL (measured at 22° C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors
(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 **TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS**

• "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

• This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

4.1 **Assay Method #1** **1,003 ± 3 µg/mL**

EDTA NIST SRM 928 Lot Number: 928

Assay Method #2 **1,002 ± 5 µg/mL**

ICP Assay NIST SRM 3159 Lot Number: 992912

13627
4/6/10

- 1.0 INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0 DESCRIPTION OF CRM** **1000 µg/mL Tin in tr. (v/v) HNO₃ / tr. (v/v) HF**
- Catalog Number: CGSN1-1, CGSN1-2, and CGSN1-5
- Lot Number: **C2-SN02024**
- Starting Material: Sn Shot
- Starting Material Purity (%): 99.991456
- Starting Material Lot No: F27N16
- Matrix: tr. (v/v) HNO₃ / tr. (v/v) HF

EXPIRES
1-May-2011**3.0 CERTIFIED VALUES AND UNCERTAINTIES****Certified Concentration:** 993 ± 4 µg/mL**Certified Density:** 0.999 g/mL (measured at 20 ± 1°C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

 (\bar{x}) = mean x_i = individual results n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

 $\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- 4.1 Assay Method #1** **993 ± 4 µg/mL (avg. of 2 runs)**
ICP Assay NIST SRM 3161a Lot Number: 070330

13638 ~~13638~~
4/6/10

1.0 **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



2.0 **DESCRIPTION OF CRM** **1000 µg/mL Titanium in 2% (v/v) HNO₃ / tr. HF**

Catalog Number: CGTI1-1, CGTI1-2, and CGTI1-5

Lot Number: **C2-TI02065**

Starting Material: Ti turnings

Starting Material Purity (%): 99.996450

Starting Material Lot No: R404TIA1

Matrix: 2% (v/v) HNO₃ / tr. HF

EXPIRES
1-May-2011

3.0 **CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 991 ± 4 µg/mL

Certified Density: 1.010 g/mL (measured at 22° C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 **TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS**

"Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

4.1 **Assay Method #1** **991 ± 4 µg/mL (avg. of 2 runs)**
 ICP Assay NIST SRM 3162a · Lot Number: 060808

Assay Method #2 **1000 ± 5 µg/mL**
 Calculated NIST SRM Lot Number: See Sec. 4.2

13030
4/6/10

- 1.0 INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0 DESCRIPTION OF CRM** **1000 µg/mL Uranium in 1.4% (v/v) HNO₃**

Catalog Number: CGU1-1, CGU1-2, and CGU1-5

Lot Number: **C2-U01079**

Starting Material: UO₂(NO₃)2.6H₂O

Starting Material Purity (%): 99.999942

Starting Material Lot No: N55947

Matrix: 1.4% (v/v) HNO₃

EXPIRES
1-May-2011

3.0 CERTIFIED VALUES AND UNCERTAINTIES

Certified Concentration: 1003 ± 3 µg/mL

Certified Density: 1.006 g/mL (measured at 20 ± 1°C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \sum_{i=1}^n x_i / n$$

\bar{x} = mean
 x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i^2)^{1/2}]}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors
(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

Certified Abundance: The ²³⁵U in this standard is depleted. The Certified abundances in Atom % are as follows:

IV's Certified Abundance	
Isotope	Atom%
Uranium 238U	99.7 ± 0.1
235U	0.29 ± 0.05

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

"Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- 4.1 Assay Method #1** **1003 ± 3 µg/mL (avg. of 2 runs)**

ICP Assay NIST SRM 3164 Lot Number: 891509

- Assay Method #2** **1000 ± 5 µg/mL**

Calculated NIST SRM Lot Number: See Sec. 4.2

13631
4/6/10

- 1.0** **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0** **DESCRIPTION OF CRM** **1000 µg/mL Vanadium in 2% (v/v) HNO₃**
- Catalog Number: CGV1-1, CGV1-2, and CGV1-5
- Lot Number: **C2-V02057**
- Starting Material: V2O₅
- Starting Material Purity (%): 99.991399
- Starting Material Lot No: Stractor 46
- Matrix: 2% (v/v) HNO₃

EXPIRES
1-May-2011

3.0 **CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 996 ± 2 µg/mL

Certified Density: 1.014 g/mL (measured at 20 ± 1°C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 **TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS**

"Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- 4.1** **Assay Method #1** **997 ± 4 µg/mL**
EDTA NIST SRM 928 Lot Number: 928
- Assay Method #2** **996 ± 2 µg/mL**
ICP Assay NIST SRM 3165 Lot Number: 992706

- 1.0 **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0 **DESCRIPTION OF CRM** **1000 µg/mL Zinc in 2% (v/v) HNO₃**
- Catalog Number: CGZN1-1, CGZN1-2, and CGZN1-5
- Lot Number: **D2-ZN02061**
- Starting Material: Zn shot
- Starting Material Purity (%): 99.999384
- Starting Material Lot No: R1207ZNA1
- Matrix: 2% (v/v) HNO₃

EXPIRES**1-May-2011**

3.0 **CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 1,007 ± 3 µg/mL

Certified Density: 1.012 g/mL (measured at 20 ± 1°C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume and the fixed error reported on the NIST SRM certificate of analysis)

4.0 **TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS**

- "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

- This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- 4.1 **Assay Method #1** **1,012 ± 3 µg/mL**
ICP Assay NIST SRM 3168a Lot Number: 080123
- Assay Method #2** **1,007 ± 3 µg/mL**
EDTA NIST SRM 928 Lot Number: 928

13595
4/6/10

tel: 800.669.6799 · 540.585.3030

fax: 540.585.3012

info@inorganicventures.com

- 1.0 INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0 DESCRIPTION OF CRM** **10000 µg/mL Aluminum in 7% v/v HNO₃**

Catalog Number: CGAL10-1, CGAL10-2, and CGAL10-5

Lot Number: **D2-AL04081**

Starting Material: Al ingot

Starting Material Purity (%): 99.998843

Starting Material Lot No: C14S012

Matrix: 7% v/v HNO₃

EXPIRES
1-May-2011

- 3.0 CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 10,015 ± 26 µg/mL

Certified Density: 1.085 g/mL (measured at 20 ± 1°C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of $k = 2$.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

- 4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS**

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- 4.1 Assay Method #1** **9996 ± 4 µg/mL**

ICP Assay NIST SRM 3101a Lot Number: 060502

- Assay Method #2** **10,015 ± 26 µg/mL**

EDTA NIST SRM 928 Lot Number: 928

1300604
4/6/10

- 1.0 INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0 DESCRIPTION OF CRM** **10000 µg/mL Calcium in 2% (v/v) HNO₃**
- Catalog Number: CGCA10-1, CGCA10-2, and CGCA10-5
- Lot Number: **C2-CA03132**
- Starting Material: CaO
- Starting Material Purity (%): 99.997483
- Starting Material Lot No: G47339
- Matrix: 2% (v/v) HNO₃

EXPIRES
1-May-2011

3.0 CERTIFIED VALUES AND UNCERTAINTIES

Certified Concentration: 10,027 ± 31 µg/mL

Certified Density: 1.040 g/mL (measured at 22° C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- 4.1 Assay Method #1** **10,027 ± 31 µg/mL**
ICP Assay NIST SRM 3109a Lot Number: 050825
- Assay Method #2** **10,048 ± 30 µg/mL**
EDTA NIST SRM 928 Lot Number: 928



13610
4/6/10

info@inorganicventures.com

- 
- Reference Material Producer
Certificate 000-00

- EXPIRES**
1-May-2011

4.1 **Assay Method #1** **9994 ± 26 µg/mL**
ICP Assay NIST SRM 3126a Lot Number: 051031

Assay Method #2 **10,011 ± 20 µg/mL**
EDTA NIST SRM 928 Lot Number: 928

13613
4/6/10

- 1.0 **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0 **DESCRIPTION OF CRM** **10000 µg/mL Magnesium in 2% (v/v) HNO₃**
- Catalog Number: CGMG10-1, CGMG10-2, and CGMG10-5
- Lot Number: **C2-MG03078**
- Starting Material: Mg metal
- Starting Material Purity (%): 99.999013
- Starting Material Lot No: 00265067003
- Matrix: 2% (v/v) HNO₃



3.0 **CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 10,038 ± 30 µg/mL

Certified Density: 1.053 g/mL (measured at 22° C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of $k = 2$.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 **TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS**

• "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

• This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term "in-house std." is specified.

- 4.1 **Assay Method #1** **10,007 ± 32 µg/mL**
ICP Assay NIST SRM 3131a Lot Number: 050302
- Assay Method #2** **10,038 ± 30 µg/mL**
EDTA NIST SRM 928 Lot Number: 928

13619
4/6/10

- 1.0 INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0 DESCRIPTION OF CRM** **10000 µg/mL Potassium in 2% (v/v) HNO₃**

Catalog Number: CGK10-1, CGK10-2, and CGK10-5

Lot Number: **C2-K03005**

Starting Material: KNO₃

Starting Material Purity (%): 99.996914

Starting Material Lot No: B19P01

Matrix: 2% (v/v) HNO₃

EXPIRES
1-May-2011

- 3.0 CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 9981 ± 28 µg/mL

Certified Density: 1.024 g/mL (measured at 22° C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of $k = 2$.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

- 4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS**

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- 4.1 Assay Method #1** **10,006 ± 52 µg/mL**

ICP Assay NIST SRM 3141a Lot Number: 051220

- Assay Method #2** **9981 ± 28 µg/mL**

Gravimetric NIST SRM Lot Number: See Sec. 4.2

13625
4610

- 1.0 **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0 **DESCRIPTION OF CRM** **10000 µg/mL Sodium in 2% HNO₃ (v/v)**
Catalog Number: CGNA10-1, CGNA10-2, and CGNA10-5
Lot Number: **B2-NA03063**
Starting Material: Na₂CO₃
Starting Material Purity (%): 99.999644
Starting Material Lot No: C181157
Matrix: 2% HNO₃ (v/v)

EXPIRES
1-May-2011

3.0 **CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 9,941 ± 42 µg/mL

Certified Density: 1.034 g/mL (measured at 22° C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 **TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS**

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- 4.1 **Assay Method #1** **9,941 ± 42 µg/mL**
 ICP Assay NIST SRM 3152a Lot Number: 010728
- Assay Method #2** **9,996 ± 12 µg/mL**
 Gravimetric NIST SRM Lot Number: See Sec. 4.2

1.0 INORGANIC VENTURES is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



2.0 DESCRIPTION OF CRM Stock Second Source Solution
Catalog No.: QCP-QCS-1
Lot Number: **D2-MEB317042**
Matrix: 5% HNO₃(v/v)

EXPIRES
1-May-2011

13646
4/10

Second Source: Whenever possible, this solution was manufactured from a second set of concentrates in our manufacturing facility.

500.00 µg/mL ea:

K, P, TI,

200.00 µg/mL ea:

As, Hg, Pb,

100.00 µg/mL ea:

Al, B, Ba, Be, Ca, Cd, Ce, Co, Cr₃, Cu, Fe, Li, Mg, Mn, Na, Ni,

Se, Sr, V, Zn,

25.00 µg/mL ea:

Ag

3.0 CERTIFIED VALUES AND UNCERTAINTIES

ELEMENT	CERTIFIED VALUE	ELEMENT	CERTIFIED VALUE	ELEMENT	CERTIFIED VALUE
Aluminum, Al	99.9 ± 0.3 µg/mL	Arsenic, As	200.0 ± 0.6 µg/mL	Barium, Ba	100.0 ± 0.1 µg/mL
Beryllium, Be	100.0 ± 0.2 µg/mL	Boron, B	100.0 ± 0.4 µg/mL	Cadmium, Cd	100.0 ± 0.2 µg/mL
Calcium, Ca	100.1 ± 0.3 µg/mL	Cerium, Ce	100.0 ± 0.3 µg/mL	Chromium+3, Cr ₃	100.1 ± 0.4 µg/mL
Cobalt, Co	100.0 ± 0.2 µg/mL	Copper, Cu	99.9 ± 0.3 µg/mL	Iron, Fe	100.0 ± 0.2 µg/mL
Lead, Pb	200.0 ± 0.8 µg/mL	Lithium, Li	100.0 ± 0.3 µg/mL	Magnesium, Mg	99.6 ± 0.4 µg/mL
Manganese, Mn	100.1 ± 0.2 µg/mL	Mercury, Hg	200.0 ± 0.5 µg/mL	Nickel, Ni	99.6 ± 0.3 µg/mL
Phosphorus, P	500.0 ± 1.5 µg/mL	Potassium, K	500.0 ± 3.0 µg/mL	Selenium, Se	100.0 ± 1.0 µg/mL
Silver, Ag	25.06 ± 0.13 µg/mL	Sodium, Na	99.8 ± 0.4 µg/mL	Strontium, Sr	100.0 ± 0.3 µg/mL
Thallium, TI	500.0 ± 2.0 µg/mL	Vanadium, V	100.0 ± 0.3 µg/mL	Zinc, Zn	99.7 ± 0.3 µg/mL

Certified Density: 1.037 g/mL (measured at 20 ± 1°C)

1364872
4/6/10

- 1.0** **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0** **DESCRIPTION OF CRM** Stock Second Source Solution

Catalog No.: QCP-QCS-2

Lot Number: **D2-MEB324018**

Matrix: tr. HF, 5% HNO₃(v/v)

EXPIRES
1-May-2011

Second Source: Whenever possible, this solution was manufactured from a second set of concentrates in our manufacturing facility.

500.00 µg/mL ea:

SiO₂, Sn,

200.00 µg/mL ea:

Sb,

100.00 µg/mL ea:

Mo, Ti

3.0 CERTIFIED VALUES AND UNCERTAINTIES

ELEMENT	CERTIFIED VALUE	ELEMENT	CERTIFIED VALUE	ELEMENT	CERTIFIED VALUE
Antimony, Sb	200.0 ± 0.5 µg/mL	Molybdenum, Mo	100.1 ± 0.4 µg/mL	Silica, SiO ₂	500.0 ± 2.4 µg/mL
Tin, Sn	500.0 ± 1.4 µg/mL	Titanium, Ti	100.0 ± 0.3 µg/mL		

Certified Density: 1.024 g/mL (measured at 20 ± 1°C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of $k = 2$.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume and the fixed error reported on the NIST SRM certificate of analysis)

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

13649
4/6/10

- 1.0** **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0** **DESCRIPTION OF CRM** Stock Second Source Solution
Catalog No.: QCP-QCS-3
Lot Number: **D2-MEB324067**
Matrix: 7% HNO₃(v/v)

EXPIRES
1-May-2011

Second Source: Whenever possible, this solution was manufactured from a second set of concentrates in our manufacturing facility.

50.00 µg/mL ea:

Se,

10.00 µg/mL ea:

Ag, Al, As, Ba, Be, Ca, Cd, Co, Cr₃, Cu, Fe, K, Mg, Mn, Mo, Na,
Ni, Pb, Sb, Th, Ti, U, V, Zn

3.0 CERTIFIED VALUES AND UNCERTAINTIES

ELEMENT	CERTIFIED VALUE	ELEMENT	CERTIFIED VALUE	ELEMENT	CERTIFIED VALUE
Aluminum, Al	10.01 ± 0.02 µg/mL	Antimony, Sb	10.02 ± 0.02 µg/mL	Arsenic, As	10.00 ± 0.03 µg/mL
Barium, Ba	10.00 ± 0.01 µg/mL	Beryllium, Be	10.00 ± 0.02 µg/mL	Cadmium, Cd	10.01 ± 0.02 µg/mL
Calcium, Ca	10.01 ± 0.03 µg/mL	Chromium+3, Cr ₃	10.00 ± 0.04 µg/mL	Cobalt, Co	10.00 ± 0.04 µg/mL
Copper, Cu	10.00 ± 0.02 µg/mL	Iron, Fe	10.01 ± 0.03 µg/mL	Lead, Pb	10.00 ± 0.02 µg/mL
Magnesium, Mg	10.01 ± 0.03 µg/mL	Manganese, Mn	10.00 ± 0.05 µg/mL	Molybdenum, Mo	10.02 ± 0.04 µg/mL
Nickel, Ni	10.00 ± 0.03 µg/mL	Potassium, K	10.01 ± 0.03 µg/mL	Selenium, Se	49.99 ± 0.13 µg/mL
Silver, Ag	10.00 ± 0.05 µg/mL	Sodium, Na	10.01 ± 0.02 µg/mL	Thallium, Tl	10.00 ± 0.02 µg/mL
Thorium, Th	9.94 ± 0.03 µg/mL	Uranium, U	10.02 ± 0.04 µg/mL	Vanadium, V	10.00 ± 0.03 µg/mL
Zinc, Zn	10.00 ± 0.04 µg/mL				

Certified Density: 1.034 g/mL (measured at 20 ± 1°C)

1.0 **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



2.0 **DESCRIPTION OF CRM** Custom Second Source Solution
Catalog No.: QCP-QCS-4
Lot Number: **C2-MEB236111**
Matrix: 7% HNO₃(v/v)

EXPIRES
1-May-2011

13654
46/40

Second Source: Whenever possible, this solution was manufactured from a second set of concentrates in our manufacturing facility.

5.00 µg/mL ea:
Hg

3.0 CERTIFIED VALUES AND UNCERTAINTIES

ELEMENT	CERTIFIED VALUE	ELEMENT	CERTIFIED VALUE	ELEMENT	CERTIFIED VALUE
Mercury, Hg	4.982 ± 0.012 µg/mL				

Certified Density: 1.034 g/mL (measured at 20 ± 1°C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

• "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

• This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

4.1 ASSAY INFORMATION

ELEMENT	METHOD	NIST SRM#	SRM LOT#	ELEMENT	METHOD	NIST SRM#	SRM LOT#
Hg	ICP Assay	3133	061204	Hg	EDTA	928	928



1.0 **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



2.0 DESCRIPTION OF CRM Stock Solution
Catalog No.: WW-LFS-1
Lot Number: **B2-MEB236093**
Matrix: 5% HNO₃(v/v)

13662
7/6/10

1,000.00 µg/mL ea:

K,

600.00 µg/mL ea:

P,

300.00 µg/mL ea:

Fe, Na,

200.00 µg/mL ea:

Al, Ce, Mg, Se, Ti,

100.00 µg/mL ea:

Ca, Pb,

80.00 µg/mL ea:

As,

70.00 µg/mL ea:

Hg,

50.00 µg/mL ea:

Ni,

40.00 µg/mL ea:

Cr₃,

30.00 µg/mL ea:

B, Cu, V,

20.00 µg/mL ea:

Ba, Be, Cd, Co, Li, Mn, Sr, Zn,

7.50 µg/mL ea:

Ag

EXPIRES

1-May-2011

3.0 CERTIFIED VALUES AND UNCERTAINTIES

ELEMENT	CERTIFIED VALUE	ELEMENT	CERTIFIED VALUE	ELEMENT	CERTIFIED VALUE
Aluminum, Al	199.5 ± 0.5 µg/mL	Arsenic, As	80.0 ± 0.3 µg/mL	Barium, Ba	20.04 ± 0.08 µg/mL
Beryllium, Be	19.96 ± 0.04 µg/mL	Boron, B	30.08 ± 0.11 µg/mL	Cadmium, Cd	20.02 ± 0.06 µg/mL
Calcium, Ca	99.9 ± 0.3 µg/mL	Cerium, Ce	200.0 ± 0.6 µg/mL	Chromium+3, Cr ₃	39.85 ± 0.15 µg/mL
Cobalt, Co	20.05 ± 0.08 µg/mL	Copper, Cu	29.92 ± 0.11 µg/mL	Iron, Fe	298.6 ± 0.6 µg/mL
Lead, Pb	100.0 ± 0.4 µg/mL	Lithium, Li	20.00 ± 0.08 µg/mL	Magnesium, Mg	199.2 ± 0.8 µg/mL
Manganese, Mn	19.90 ± 0.05 µg/mL	Mercury, Hg	69.9 ± 0.1 µg/mL	Nickel, Ni	49.87 ± 0.14 µg/mL

1.0 INORGANIC VENTURES is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



2.0 DESCRIPTION OF CRM Stock Solution
Catalog No.: WW-LFS-2
Lot Number: **B2-MEB236102**
Matrix: tr. HF, 5% HNO₃(v/v)

13661
4/6/10

200.00 µg/mL ea:
SiO₂,
80.00 µg/mL ea:
Sb,
70.00 µg/mL ea:
Sn,
40.00 µg/mL ea:
Mo,
20.00 µg/mL ea:
Ti

EXPIRES
1-May-2011

3.0 CERTIFIED VALUES AND UNCERTAINTIES

ELEMENT	CERTIFIED VALUE	ELEMENT	CERTIFIED VALUE	ELEMENT	CERTIFIED VALUE
Antimony, Sb	80.2 ± 0.2 µg/mL	Molybdenum, Mo	40.05 ± 0.10 µg/mL	Silica, SiO ₂	199.5 ± 0.4 µg/mL
Tin, Sn	69.9 ± 0.2 µg/mL	Titanium, Ti	20.04 ± 0.04 µg/mL		

Certified Density: 1.022 g/mL (measured at 22° C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of $k = 2$.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i^2)^{1/2}]}{(n)^{1/2}}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

1.0 **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



13643
4/6/10

2.0 DESCRIPTION OF CRM Stock Solution
Catalog No.: 6020ICS-0A
Lot Number: **C2-MEB307012**
Matrix: 1.4% HNO₃(v/v)

EXPIRES
1-May-2011

10,000.00 µg/mL ea:

Chloride,

2,000.00 µg/mL ea:

C,

1,000.00 µg/mL ea:

Al, Ca, Fe, K, Mg, Na, P, S,

20.00 µg/mL ea:

Mo, Ti

3.0 CERTIFIED VALUES AND UNCERTAINTIES

ELEMENT	CERTIFIED VALUE	ELEMENT	CERTIFIED VALUE	ELEMENT	CERTIFIED VALUE
Aluminum, Al	1,002 ± 3 µg/mL	Calcium, Ca	1,002 ± 3 µg/mL	Carbon, C	2,004 ± 10 µg/mL
Chloride, Chloride	10,020.0 ± 20.0 µg/mL	Iron, Fe	1,002 ± 3 µg/mL	Magnesium, Mg	1,002 ± 3 µg/mL
Molybdenum, Mo	20.04 ± 0.07 µg/mL	Phosphorus, P	1,002 ± 5 µg/mL	Potassium, K	1,002 ± 3 µg/mL
Sodium, Na	1,002 ± 2 µg/mL	Sulfur, S	1,002 ± 5 µg/mL	Titanium, Ti	20.09 ± 0.06 µg/mL

Certified Density: 1.033 g/mL (measured at 22° C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

Certified Value (\bar{x}) = $\frac{\sum x_i}{n}$

Uncertainty (\pm) = $\frac{2(\sum s_i^2)^{1/2}}{(n)^{1/2}}$

(\bar{x}) = mean


x_i = individual results

n = number of measurements

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

1.0 INORGANIC VENTURES is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



2.0 DESCRIPTION OF CRM Stock Solution
Catalog No.: 6020ICS-0B
Lot Number: **A2-MEB194156**
Matrix: 2% HNO₃(v/v)

EXPIRES
1-May-2011

13645
4/6/10

2.00 µg/mL ea:

Ag, As, Cd, Co, Cr₃, Cu, Mn, Ni, Zn

3.0 CERTIFIED VALUES AND UNCERTAINTIES

ELEMENT	CERTIFIED VALUE	ELEMENT	CERTIFIED VALUE	ELEMENT	CERTIFIED VALUE
Arsenic, As	2.002 ± 0.012 µg/mL	Cadmium, Cd	2.004 ± 0.008 µg/mL	Chromium+3, Cr ₃	2.001 ± 0.004 µg/mL
Cobalt, Co	2.000 ± 0.008 µg/mL	Copper, Cu	1.998 ± 0.006 µg/mL	Manganese, Mn	2.005 ± 0.004 µg/mL
Nickel, Ni	2.004 ± 0.004 µg/mL	Silver, Ag	2.004 ± 0.006 µg/mL	Zinc, Zn	1.998 ± 0.004 µg/mL

Certified Density: 1.009 g/mL (measured at 22° C)

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2[(\sum S_i)^2]^{1/2}}{(n)^{1/2}}$$

(\bar{x}) = mean
 x_i = individual results
 n = number of measurements
 $\sum S_i$ = The summation of all significant estimated errors
(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.



1.0 INORGANIC VENTURES is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



2.0 DESCRIPTION OF CRM Stock Solution
Catalog No.: 6020SPK-W
Lot Number: **B2-CICP22006**
Matrix: 5% HNO₃(abs)

13644
4/6/10

EXPIRES
1-May-2011

100.00 µg/mL ea:

Fe,

50.00 µg/mL ea:

Ba, Zn,

20.00 µg/mL ea:

Co, Cr₃, Cu, Mn, Ni, Sb, V,

10.00 µg/mL ea:

As, Pb,

5.00 µg/mL ea:

Ag, Be, Cd, Se, Ti

3.0 CERTIFIED VALUES AND UNCERTAINTIES

ELEMENT	CERTIFIED VALUE	ELEMENT	CERTIFIED VALUE	ELEMENT	CERTIFIED VALUE
Antimony, Sb	20.03 ± 0.04 µg/mL	Arsenic, As	10.04 ± 0.03 µg/mL	Barium, Ba	50.15 ± 0.08 µg/mL
Beryllium, Be	4.989 ± 0.010 µg/mL	Cadmium, Cd	5.016 ± 0.015 µg/mL	Chromium+3, Cr ₃	19.99 ± 0.08 µg/mL
Cobalt, Co	20.01 ± 0.10 µg/mL	Copper, Cu	20.04 ± 0.04 µg/mL	Iron, Fe	99.7 ± 0.2 µg/mL
Lead, Pb	10.04 ± 0.02 µg/mL	Manganese, Mn	20.01 ± 0.04 µg/mL	Nickel, Ni	20.02 ± 0.06 µg/mL
Selenium, Se	5.001 ± 0.020 µg/mL	Silver, Ag	5.015 ± 0.015 µg/mL	Thallium, Tl	5.013 ± 0.010 µg/mL
Vanadium, V	20.01 ± 0.04 µg/mL	Zinc, Zn	50.00 ± 0.18 µg/mL		

Certified Density: 1.040 g/mL (measured at 22° C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

13642
4/6/10

1.0 INORGANIC VENTURES is an ISO Guide 34:2000 registered Certified Reference Material (CRM) Manufacturer (Certificate #883-02). The certificate is designed and the data is determined in accordance with ISO Guide 31:2000 (Reference Materials-Contents of Certificates and Labels), ISO Guide 34:2000 "Quality System Guidelines for the Production of Reference Materials," and ISO Guide 35-1989 "Certification of Reference Materials - General and Statistical Principles."

2.0 DESCRIPTION OF CRM Stock Solution
Catalog No.: 6020SPK-S
Lot Number: **A2-CICP20089**
Matrix: 5% HNO₃(abs)

50.00 µg/mL each:
Ba, Cr₃, Cu, Zn,
30.00 µg/mL each:
V,
25.00 µg/mL each:
Ni,
20.00 µg/mL each:
Co, Pb, Sb,
10.00 µg/mL each:
Ag, As, Cd,
5.00 µg/mL each:
Be, Se, Ti

EXPIRES
1-May-2011

3.0 CERTIFIED VALUES AND UNCERTAINTIES

ELEMENT	CERTIFIED VALUE	ELEMENT	CERTIFIED VALUE	ELEMENT	CERTIFIED VALUE
Antimony, Sb	20.07 ± 0.06 µg/mL	Arsenic, As	10.01 ± 0.03 µg/mL	Barium, Ba	50.10 ± 0.09 µg/mL
Beryllium, Be	5.007 ± 0.020 µg/mL	Cadmium, Cd	10.02 ± 0.02 µg/mL	Chromium+3, Cr ₃	50.11 ± 0.06 µg/mL
Cobalt, Co	20.03 ± 0.09 µg/mL	Copper, Cu	49.88 ± 0.16 µg/mL	Lead, Pb	20.04 ± 0.06 µg/mL
Nickel, Ni	25.00 ± 0.09 µg/mL	Selenium, Se	5.020 ± 0.010 µg/mL	Silver, Ag	10.00 ± 0.03 µg/mL
Thallium, Tl	5.010 ± 0.015 µg/mL	Vanadium, V	29.92 ± 0.09 µg/mL	Zinc, Zn	49.93 ± 0.13 µg/mL

Certified Density: 1.036 g/mL (measured at 22° C)

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty:

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

1.0 INORGANIC VENTURES is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



13639
4/6/10

2.0 DESCRIPTION OF CRM Stock Solution
Catalog No.: 2007ICS-1
Lot Number: **A2-MEB236039**
Matrix: tr. HF, 2% HNO₃(abs)

1,000.00 µg/mL ea:
Ti,
500.00 µg/mL ea:
B,
300.00 µg/mL ea:
Mo,
230.00 µg/mL ea:
Si

EXPIRES
1-May-2011

3.0 CERTIFIED VALUES AND UNCERTAINTIES

ELEMENT	CERTIFIED VALUE	ELEMENT	CERTIFIED VALUE	ELEMENT	CERTIFIED VALUE
Boron, B	499.9 ± 2.1 µg/mL	Molybdenum, Mo	299.9 ± 0.5 µg/mL	Silicon, Si	230.0 ± 0.5 µg/mL
Titanium, Ti	1,000 ± 2 µg/mL				

Certified Density: 1.019 g/mL (measured at 22° C)

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2\sqrt{\frac{\sum (x_i - \bar{x})^2}{n-1}}}{n^{1/2}}$$

(\bar{x}) = mean
 x_i = individual results
n = number of measurements
 $\sum x_i$ = The summation of all significant estimated errors
(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

- "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)
- This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

13640
4/6/10

- 1.0 **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0 **DESCRIPTION OF CRM** Custom Solution
Catalog No.: 2007ICS-3
Lot Number: **C2-MEB307125**
Matrix: 7% HNO₃(v/v)

EXPIRES
1-May-2011

20,000.00 µg/mL ea:

K,

1,000.00 µg/mL ea:

As, Pb, Tl,

500.00 µg/mL ea:

Se,

300.00 µg/mL ea:

Ag, Ba, Cd, Co, Cr₃, Cu, Ni, V, Zn,

200.00 µg/mL ea:

Mn,

100.00 µg/mL ea:


Be

3.0 CERTIFIED VALUES AND UNCERTAINTIES

ELEMENT	CERTIFIED VALUE	ELEMENT	CERTIFIED VALUE	ELEMENT	CERTIFIED VALUE
Arsenic, As	1,000 ± 3 µg/mL	Barium, Ba	300.0 ± 0.6 µg/mL	Beryllium, Be	100.2 ± 0.4 µg/mL
Cadmium, Cd	300.0 ± 0.7 µg/mL	Chromium+3, Cr ₃	300.0 ± 1.1 µg/mL	Cobalt, Co	300.0 ± 1.2 µg/mL
Copper, Cu	300.8 ± 0.8 µg/mL	Lead, Pb	1,000 ± 2 µg/mL	Manganese, Mn	199.3 ± 0.4 µg/mL
Nickel, Ni	300.2 ± 0.8 µg/mL	Potassium, K	20,000.0 ± 50.0 µg/mL	Selenium, Se	500.0 ± 2.0 µg/mL
Silver, Ag	300.0 ± 0.7 µg/mL	Thallium, Tl	1,000 ± 2 µg/mL	Vanadium, V	300.0 ± 0.9 µg/mL
Zinc, Zn	300.2 ± 1.3 µg/mL				

Certified Density: 1.089 g/mL (measured at 22° C)

1.0 INORGANIC VENTURES is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



2.0 DESCRIPTION OF CRM Stock Solution
Catalog No.: 2007ICS-4
Lot Number: **B2-MEB236099**
Matrix: 3% HNO₃(v/v)

EXPIRES
1-May-2011

15,000.00 µg/mL ea:
Ca,
12,500.00 µg/mL ea:
Fe,
7,500.00 µg/mL ea:
Mg,
3,000.00 µg/mL ea:
Al,
2,500.00 µg/mL ea:
Na

3.0 CERTIFIED VALUES AND UNCERTAINTIES

ELEMENT	CERTIFIED VALUE	ELEMENT	CERTIFIED VALUE	ELEMENT	CERTIFIED VALUE
Aluminum, Al	3,001 ± 8 µg/mL	Calcium, Ca	15,000.0 ± 40.0 µg/mL	Iron, Fe	12,500.0 ± 20.0 µg/mL
Magnesium, Mg	7,500.0 ± 30.0 µg/mL	Sodium, Na	2,502 ± 11 µg/mL		

Certified Density: 1.180 g/mL (measured at 22° C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

(\bar{x}) = mean
 x_i = individual results
n = number of measurements
 $\sum s_i$ = The summation of all significant estimated errors
(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

4.0 TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

- "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)
- This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.



TechLaw, Inc.
Environmental Services Assistance Team
16194 W. 45th Drive, Golden, CO 80403
303-312-7720

Task Order: 32
Contract: EP-W-06-33
Valid: April 2010 – May 2011

Certificates of Analysis

Valid through May 1, 2010

Alkalinity

Anions by Ion Chromatography

Dissolved Organic Carbon

- Initial Calibration Verification (ICV) Standards
- Laboratory Check Standards (LCS)
- Matrix Spike Solutions
- Interference Check (ICSA / AB) Standards

- 1.0 **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0 **DESCRIPTION OF CRM** Ion Chromatography 10,000 µg/mL Sulfate in H₂O
- Catalog Number: ICSO410-1 and ICSO410-5
- Lot Number: **B2-SOX01084**
- Starting Material: K₂SO₄
- Starting Material Purity (%): 99.0000
- Starting Material Lot No.: 08524KC
- Matrix: H₂O

13657
4/6/10**EXPIRES**
1-May-2011

3.0 **CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 10,030 ± 22 µg/mL

Certified Density: 1.011 g/mL (measured at 20 ± 1°C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of $k = 2$.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume and the fixed error reported on the NIST SRM certificate of analysis)

4.0 **TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS**

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

Assay Method #1 10,030 ± 22 µg/mL (Avg. of 2 Runs)
IC Assay NIST SRM 3154 Lot Number: 892205

Assay Method #2 10,010 ± 50 µg/mL
Calculated NIST SRM Lot Number: See Sec. 4.2

CERTIFICATE OF ANALYSIS

tel: 800.669.6799 · 540.585.3030
fax: 540.585.3012
info@inorganicventures.com

1.0

INORGANIC VENTURES is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



2.0

DESCRIPTION OF CRM

Ion Chromatography 10000 µg/mL Nitrate as N in H₂O
ICNNO310-1 and ICNNO310-5

Catalog Number:

B2-NOX02064

Lot Number:

NaNO₃

Starting Material:

Starting Material Purity (%): **99.0000**

Starting Material Lot No.: **12616AC**

Matrix:

H₂O

13660
4660

3.0

CERTIFIED VALUES AND UNCERTAINTIES

Certified Concentration: 10,056 ± 30 µg/mL

Certified Density:

1.038 g/mL (measured at 20 ± 1°C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of $k = 2$.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume and the fixed error reported on the NIST SRM certificate of analysis)

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS

"Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

Assay Method #1

10,056 ± 30 µg/mL (Avg. of 2 Runs)
IC Assay NIST SRM 3185 Lot Number: 050517

Assay Method #2

10,008 ± 50 µg/mL

Calculated NIST SRM Lot Number: See Sec. 4.2

EXPIRES
1-May-2011

- 1.0** **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0** **DESCRIPTION OF CRM** **Ion Chromatography 10000 µg/mL Nitrite as N in H₂O**
- Catalog Number: ICNNO210-1 and ICNNO210-5
- Lot Number: **B2-NOX02059**
- Starting Material: NaNO₂
- Starting Material Purity (%): 99.6000
- Starting Material Lot No.: 18122HO
- Matrix: H₂O

13658
4/6/10

EXPIRES
1-May-2011

3.0 **CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 9,984 ± 30 µg/mL

Certified Density: 1.037 g/mL (measured at 20 ± 1°C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume and the fixed error reported on the NIST SRM certificate of analysis)

4.0 **TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS**

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

Assay Method #1 **9,984 ± 30 µg/mL (Avg. of 2 Runs)**
 IC Assay NIST SRM 40h Lot Number: 000412

Assay Method #2 **10,009 ± 50 µg/mL**
 Calculated NIST SRM Lot Number: See Sec. 4.2

- 1.0 **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0 **DESCRIPTION OF CRM** Ion Chromatography 10,000 µg/mL Chloride in H₂O
- Catalog Number: ICCL10-1 and ICCL10-5
- Lot Number: **B2-CL01086**
- Starting Material: KCl
- Starting Material Purity (%): 99.9990
- Starting Material Lot No.: 075K0024
- Matrix: H₂O

13659
4/6/10**EXPIRES**
1-May-2011

3.0 **CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 10,202 ± 20 µg/mL

Certified Density: 1.010 g/mL (measured at 20 ± 1°C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of $k = 2$.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume and the fixed error reported on the NIST SRM certificate of analysis)

4.0 **TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS**

· "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

· This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- Assay Method #1** **10,202 ± 20 µg/mL**
IC Assay NIST SRM 3182 Lot Number: 060925
- Assay Method #2** **10,180 ± 40 µg/mL**
Volhard NIST SRM 999b Lot Number: 999b

13601
4/6/10

tel: 800.669.6799 · 540.585.3030

fax: 540.585.3012

info@inorganicventures.com

- 1.0** **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0** **DESCRIPTION OF CRM** **10000 µg/mL Bromide (NH₄) in H₂O**
- Catalog Number: CGICBR10-1 and CGICBR10-5
- Lot Number: **A2-BR01066**
- Starting Material: NH₄Br
- Starting Material Purity (%): 99.998491
- Starting Material Lot No: DI05205EU
- Matrix: H₂O

Exp 01-MAY-2011

3.0 **CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 9993 ± 34 µg/mL

Certified Density: 1.004 g/mL (measured at 20 ± 1°C)

The Certified Value is based upon the most precise method used to analyze this CRM. The following equations are used in the calculation of the certified value and the uncertainty.

Certified Value $(\bar{x}) = \frac{\sum x_i}{n}$ (\bar{x}) = mean

x_i = individual results

n = number of measurements

Uncertainty $(\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume and the fixed error reported on the NIST SRM certificate of analysis)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of $k = 2$.

4.0 **TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS**

"Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

- 4.1** **Assay Method #1** **9993 ± 34 µg/mL**
IC Assay NIST SRM 3184 Lot Number: 020701
- Assay Method #2** **9976 ± 47 µg/mL**
Volhard NIST SRM 999b Lot Number: 999b

- 1.0 **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0 **DESCRIPTION OF CRM** Ion Chromatography 10,000 µg/mL Fluoride in H₂O
- Catalog Number: ICF10-1 and ICF10-5
- Lot Number: **C2-F01054**
- Starting Material: Sodium Fluoride
- Starting Material Purity (%): 99.0000
- Starting Material Lot No.: 04408EZ
- Matrix: H₂O

13667
7/6/10**EXPIRES**
1-May-2011

3.0 **CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 10,117 ± 19 µg/mL

Certified Density: 1.020 g/mL (measured at 20 ± 1°C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume and the fixed error reported on the NIST SRM certificate of analysis)

4.0 **TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS**

• "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

• This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

Assay Method #1 10,117 ± 19 µg/mL (Avg. of 2 Runs)
IC Assay NIST SRM 3183 Lot Number: 050721

Assay Method #2 10,019 ± 51 µg/mL
Calculated NIST SRM Lot Number: See Sec. 4.2

- 1.0 **INORGANIC VENTURES** is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



- 2.0 **DESCRIPTION OF CRM** Ion Chromatography 10,000 µg/mL Phosphate as P in H₂O
- Catalog Number: ICPPO410-1 and ICPPO410-5
- Lot Number: **B2-POX01082**
- Starting Material: NH₄H₂PO₄
- Starting Material Purity (%): 99.9990
- Starting Material Lot No.: 10911CA / 07430AE
- Matrix: H₂O

13665
4/6/10**EXPIRES**
1-May-2011

3.0 **CERTIFIED VALUES AND UNCERTAINTIES**

Certified Concentration: 10,006 ± 26 µg/mL

Certified Density: 1.018 g/mL (measured at 20 ± 1°C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume and the fixed error reported on the NIST SRM certificate of analysis)

4.0 **TRACEABILITY TO NIST AND VALUES OBTAINED BY INDEPENDENT METHODS**

• "Property of the result of a measurement or the value of a standard whereby it can be related to stated references, usually national or international standards, through an unbroken chain of comparisons all having stated uncertainties." (ISO VIM, 2nd ed., 1993, definition 6.10)

• This product is Traceable to NIST via an unbroken chain of comparisons. The uncertainties for each certified value are reported, taking into account the SRM uncertainty error and the measurement, weighing and volume dilution errors. In rare cases where no NIST SRMs are available, the term 'in-house std.' is specified.

Assay Method #1 **10,006 ± 26 µg/mL (Avg. of 2 Runs)**
IC Assay NIST SRM 3186 Lot Number: 000330

Assay Method #2 **10,009 ± 50 µg/mL**
Calculated NIST SRM Lot Number: See Sec. 4.2

1.0 INORGANIC VENTURES is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



2.0 DESCRIPTION OF CRM Stock Second Source Custom Solution
Catalog No.: QCP-QCS-5
Lot Number: **C2-MEB307127**
Matrix: H2O

13653
46/10

Second Source: Whenever possible, this solution was manufactured from a second set of concentrates in our manufacturing facility.

75.00 µg/mL ea:

Sulfate,

50.00 µg/mL ea:

Bromide,

25.00 µg/mL ea:

oPhosph

ate_as_

P,

15.00 µg/mL ea:

Chloride, Nitrite_a

s_N,

10.00 µg/mL ea:

Fluoride, Nitrate_a

s_N

EXPIRES
1-May-2011

3.0 CERTIFIED VALUES AND UNCERTAINTIES

ELEMENT	CERTIFIED VALUE	ELEMENT	CERTIFIED VALUE	ELEMENT	CERTIFIED VALUE
Bromide, Bromide	49.98 ± 0.24 µg/mL	Chloride, Chloride	15.00 ± 0.04 µg/mL	Fluoride, Fluoride	10.00 ± 0.03 µg/mL
Nitrate_as_N, Nitrate_as_N	10.02 ± 0.03 µg/mL	Nitrite_as_N, Nitrite_as_N	15.00 ± 0.04 µg/mL	o-Phosphate as P, oPhosphate_25.00	± 0.07 µg/mL
Sulfate, Sulfate	75.0 ± 0.2 µg/mL				

Certified Density: 0.997 g/mL (measured at 20 ± 1°C)

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean

x_i = individual results

n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2f(\sum s_i)^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

INORGANIC VENTURES is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."


1365/
4/6/10

WATER QC CERTIFIED REFERENCE MATERIAL Complex Nutrients

Catalog No: QCP-NUT-2

Lot Number: B2-NUT01111

EXPIRES

1-May-2011

STABILITY AND STORAGE INFORMATION - This CRM can be stored at room temperature before opening. After opening and dilution, the EPA recommends that it be stored at 4 °C for no more than 48 hours for Nitrate as N and Phosphate as P. For Nitrate plus Nitrite as N and Ammonia as N samples, a "maximum" holding time of 28 days at 4 °C is recommended. The EPA recommendations for holding time and storage conditions should be followed after opening and dilution.

SPECIFICATIONS AND TRACEABILITY:

Parameter	Certified Value	Made to Value	Analytical Method	NIST Traceability	Acceptance Limits
Total Kjeldahl Nitrogen as N	15.1 ± 0.3 mg/L	15.0 mg/L	Dumas Micro Combustion	SRM 141c	18.6 – 11.5 mg/L
Total Organic Phosphorus as P	8.02 ± 0.04 mg/L	8.00 mg/L	ICP	SRM 3139a	9.64 – 6.40 mg/L

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

 (\bar{x}) = mean

 x_i = individual results

 n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

 $\sum s_i$ = The summation of all significant estimated errors

(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis)

ANALYZED DENSITY OF SOLUTION (measured at 22°C): 1.013 g/mL

INORGANIC VENTURES is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



WATER QC Reference Material Solids
Catalog No: QCP-SLD

EXPIRES
1-May-2011

Lot Number: **C2-SLD02007**

13647
4/6/10

STABILITY AND STORAGE INFORMATION - This reference material can be stored at room temperature before and after opening. The EPA recommends a "maximum" holding time for solids samples of 7 days at 4 °C. Our stability data indicates that this standard should be disposed of in 3 months after opening.

SPECIFICATIONS AND TRACEABILITY:

Parameter	Certified Value	Made to Value	Analytical Method	NIST Traceability	Acceptance Limits
Filterable Residue	3893.3 ± 68 mg/L	4200 mg/L	EPA Method 160.1	Gravimetric	4835 - 2952 mg/L
Non-filterable Residue	151.6 ± 3.2 mg/L	185 mg/L	EPA Method 160.2	Gravimetric	161.9 - 141.3 mg/L
Total Residue	4052 ± 116 mg/L	4385 mg/L	EPA Method 160.3	Gravimetric	5032 - 3073 mg/L

The following equations are used in the calculation of the certified value and the uncertainty. Reported uncertainties represent expanded uncertainties expressed at approximately the 95% confidence level using a coverage factor of k = 2.

$$\text{Certified Value } (\bar{x}) = \frac{\sum x_i}{n}$$

(\bar{x}) = mean
 x_i = individual results
n = number of measurements

$$\text{Uncertainty } (\pm) = \frac{2[(\sum s_i)^2]^{1/2}}{(n)^{1/2}}$$

$\sum s_i$ = The summation of all significant estimated errors
(Most common are the errors from instrumental measurement, weighing, dilution to volume, and the fixed error reported on the NIST SRM certificate of analysis.)

INORGANIC VENTURES is an ISO Guide 34 "General Requirements for the Competence of Reference Material Producers" and ISO 9001:2000 registered manufacturer. Our manufacturing laboratory is accredited to ISO/IEC 17025 "General Requirements for the Competence of Testing and Calibration Laboratories."



Lims:
0101501

WATER QC REFERENCE MATERIAL Minerals

Catalog No: QCP-MIN

EXPIRES
01/2011

Lot Number: C2-MIN01121

STABILITY AND STORAGE INFORMATION

- Do not put transfer devices, probes, etc. in sample container. The insertion of a pH electrode, for example, can significantly increase the conductivity, potassium, and chloride values. This standard can be stored at room temperature before opening. After opening, the EPA recommends a "maximum" holding time for the following:

PARAMETER	HOLDING CONDITIONS	HOLDING TIME
Alkalinity	4°C	14 days
Conductivity	4°C	28 days
Chloride	None required	28 days
Sulfate	4°C	28 days
Nitrate as N	4°C	48 hours
Fluoride	None required	28 days
Sodium	HNO ₃ to pH<2	6 months
Potassium	HNO ₃ to pH<2	6 months

*pH: The value listed below is for informational purposes only. The pH value of this CRM is not stable and cannot be relied upon. It can change up to 1 pH unit. For a certified pH CRM, use catalog no. QCP-PH.

SPECIFICATIONS AND TRACEABILITY:

Parameter	Certified Value	Made to Value	Analytical Method	NIST Traceability	Acceptance Limits
Alkalinity	123.01 ± 0.73 mg/L CaCO ₃	123.7 mg/L	EPA Method 310.1	723d	130.84 – 115.18 mg/L CaCO ₃
Conductivity	1186 ± 1 µmhos/cm @ 25°C	Measured	EPA Method 120.1	999b	1310 – 1062 µmhos/cm @ 25°C
Chloride	192.211 ± 1.601 mg/L	198.0 mg/L	EPA Method 300.0	3182	206.496 – 177.926 mg/L
Fluoride	5.633 ± 0.177 mg/L	6.000 mg/L	EPA Method 300.0	3183	6.089 – 5.176 mg/L
Sulfate	113.542 ± 2.974 mg/L	120.0 mg/L	EPA Method 300.0	3154	129.612 – 97.472 mg/L
Nitrate as N	4.685 ± 0.194 mg/L	5.001 mg/L	EPA Method 300.0	3185	5.635 – 3.735 mg/L
Sodium	193.904 ± 6.657 mg/L	242.7 mg/L	ICP	3152a	212.83 – 174.985 mg/L
Potassium	92.095 ± 1.733 mg/L	97.63 mg/L	ICP	3141a	105.144 – 79.047 mg/L
pH	9.18 units	Measured	EPA Method 310.1	186g, 185h	*See parameters table above



RICCA CHEMICAL COMPANY

Arlington, TX 76012

Pocomoke City, MD 21851

Batesville, IN 47006

<http://www.riccachemical.com>

1-888-GO-RICCA

customerservice@riccachemical.com

Certificate of Analysis

Conductivity/TDS Standard, 1 mS/cm (1000 µmho/cm) at 25°C, 495 ppm as NaCl

Lot Number: 1909460

Product Number: 2243

Expiration Date: MAR 2011

Manufacture Date: 9/18/2009

The certified value for this product is confirmed in independent testing by a second qualified chemist.

Contains:

Name	CAS#	Grade
Sodium Chloride, NaCl	7647-14-5	ACS
Water, Deionized, H ₂ O	7732-18-5	ACS, ASTM D 1193 (Type I), EP, USP

Test Name	Assay Method	Specification	Result
Appearance	Clarity, Color, Odor	Clear, colorless, odorless	Passed Test
Conductivity at 25 °C (traceable to NIST SRM 3193)	Conductivity determination	1000 ± 1 µS/cm at 25.0 °C	1000 µS/cm at 25.0 °C

Volumetric glassware complies with Class A tolerance requirements of ASTM E 288 and NIST Circular 434; it is calibrated before first use and recalibrated regularly in accordance with ASTM E 542 and NIST Procedure NBSIR 74-481. Balances are calibrated regularly with weights certified traceable to the NIST national mass standard. Thermometers and temperature probes are calibrated before first use and recalibrated regularly with a thermometer traceable to NIST standards. All products are prepared according to master documents that assure manufacture according to validated methods. Batch records document raw material traceability and production and testing history for each lot manufactured.

Shelf Life (unopened container):

Part Number	Shelf Life
2243-32	18 months
2243-5	18 months
2243-20	18 months
2243-1	18 months
2243-16	18 months
2243-100P	18 months

Recommended Storage: 15°C - 30°C (59°F - 86°F)

LaNelle Ohlhausen
Quality Assurance

This Certificate of Analysis is designed to comply with ISO Guide 31 "Reference Materials -- Contents of Certificates and Labels."

To determine manufacturing site using lot number, visit www.riccachemical.com/AboutUs/lot.pdf.



Certificate of Analysis

Analytical Solutions

Total Organic Carbon (TOC) Standard

Catalog Number: IQC-106

Lot Number: K00587

Job Number: J00009977

Lot Issue Date: 06/17/2009

Expiration Date: 07/31/2011

This Certified Reference Material (CRM) was manufactured and verified in accordance with ULTRA's ISO 9001 registered quality system. The analyte concentrations were verified by our ISO 17025 accredited laboratory to be within $\pm 2.5\%$, when compared to calibration standards independently prepared using NIST SRM(s). The certified value and uncertainty value for each analyte is determined gravimetrically.

Analyte	True Value			Analytical Method	NIST SRM
TOC	1000	\pm	5 mg/L	TOC Analyzer	84K

Matrix: low TOC water (< 50 ppb)

ULTRA uses purified acids, 18 megohm double deionized water, calibrated Class A glassware & meticulously cleaned bottles in the manufacturing of ULTRAgrade standards. Balances used in the manufacturing of this standard are calibrated with weights traceable to NIST in compliance with ANSI/NCSL Z-540-1 and ISO 9001



ISO 17025:2005
Accredited
A2LA
Cert. No. 0851.01

ISO 9001:2000
Registered
TUV USA, Inc.
Cert. No. 06-1004

250 Smith Street, North Kingstown, RI 02852 USA

401-294-9400 Fax: 401-295-2330

www.ultrasci.com

See Reverse For Additional Information

William J. Leary
Quality Assurance Manager



Alkalinity - EPA Method 310.1

Mettler Auto Titrator

Date: 11/30/2010

Analyst: Scott Van Overmolen

Work Order(s): C101104

TDF DE-220

Standards Information

CCV LIMS #: 9072201 Concentration 100 mg/L Exp. Date 12/22/2010

ICV/LCS: QCP-MIN LIMS ID#: 0101501 Concentration 123.01 mg/L Exp. Date 11/1/2011

Solution Information

0.1 N H₂SO₄ Stock Solution Lot #: C08506

~0.0125g per 50 mL; i.e. 1mL of (1.25g / 100mL) to 50 mL for titer determinations.

Date Titer Solutions Prepared: 11/30/2010

0.02 N H₂SO₄ Titer Value: 0.02

Analysis Information

pH Calibration Slope: -58.10

Titer Mean (N): 0.01973

Number of samples: 33

QC Sample ID: C101104-01
C101104-52

Comments / Information

Sequence: 1012001/1012002

Batch: 101115/101116

LIMS Upload 12/01/2010

Peer Review By: SW Date: 12-2-10

PREPARATION BENCH SHEET

1011115

TechLaw, Inc. - ESAT Region 8

Matrix: Water

Date Prepared: 11/30/10 08:21 By: SV

Prepared using: ADMIN - No Lab Prep Req'd

Printed: 11/30/2010 8:23:03AM

Lab Number	Analysis	EPA Tag ID	Initial (mL)	Final (mL)	Spike1 ID	ul Spike1	Spike2 ID	ul Spike2	Source ID	QC Code	Extraction Comments
1011115-BLK1	QC		50	50						Blank	
1011115-DUP1	QC		50	50					C101104-01	Duplicate	
1011115-SRM1	QC		50	50	0101501	50000				Reference	
C101104-01 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLDRBG	
C101104-04 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLDRBG DUP	
C101104-07 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLDRMZ1a	
C101104-10 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLDRMZ1b	
C101104-13 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLDRMZ1c	
C101104-16 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLDRMZ2	
C101104-19 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLPO01	
C101104-22 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLPO02	
C101104-25 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLPO03	
C101104-28 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLPO04	
C101104-31 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLPO05	
C101104-34 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLSW01	
C101104-37 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLSW02	
C101104-40 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLSW03	
C101104-43 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLSW04	
C101104-46 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLSW05	
C101104-49 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLSWDR3	

DG-220

PREPARATION BENCH SHEET

1011116

TechLaw, Inc. - ESAT Region 8

Printed: 11/30/2010 8:23:44AM

Matrix: Water

Date Prepared: 11/30/10 08:23 By: SV

Prepared using: ADMIN - No Lab Prep Req'd

Lab Number	Analysis	EPA Tag ID	Initial (mL)	Final (mL)	Spike1 ID	ul Spike1	Spike2 ID	ul Spike2	Source ID	QC Code	Extraction Comments
1011116-BLK1	QC		50	50						Blank	
1011116-DUP1	QC		50	50					C101104-52	Duplicate	
1011116-SRM1	QC		50	50	0101501	50000				Reference	
C101104-52 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLSWDR4	
C101104-55 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLSWDR6	
C101104-58 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLSWDR7b	
C101104-61 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLSWDR7c	
C101104-64 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLSWFB	
C101104-67 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLSWP06	
C101104-70 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLSWP07a	
C101104-73 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLSWP07b	
C101104-76 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLSWP08	
C101104-79 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLSWP09	
C101104-82 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLSWP10	
C101104-85 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLSWP11	
C101104-88 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLSWP12	
C101104-91 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLSWP14	
C101104-94 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLSWP15	
C101104-97 A	WC - Alkalinity	lo Tag Prefix-	50	50						SLSWPP	

DG-220

ANALYSIS SEQUENCE

1012001

11/30/10

Instrument: Mettler AT

Sequence Date: 11/30/10 00:00

Printed: 12/1/2010 9:54:15AM

Lab Number	Dilut. Factor	Analysis	STD ID	Sample/Std Name	EPA Tag ID	Source Sple	Comments
1011115-SRM1		QC		Reference		-	
1011115-BLK1		QC		Blank		-	
C101104-01 A		WC - Alkalinity		SLDRBG	No Tag Prefix-A		
1011115-DUP1		QC		Duplicate		C101104-01	
C101104-04 A		WC - Alkalinity		SLDRBG DUP	No Tag Prefix-A		
C101104-07 A		WC - Alkalinity		SLDRMZ1a	No Tag Prefix-A		
C101104-10 A		WC - Alkalinity		SLDRMZ1b	No Tag Prefix-A		
C101104-13 A		WC - Alkalinity		SLDRMZ1c	No Tag Prefix-A		
C101104-16 A		WC - Alkalinity		SLDRMZ2	No Tag Prefix-A		
C101104-19 A		WC - Alkalinity		SLPO01	No Tag Prefix-A		
1012001-CCV1		QC	9072201	Calibration Check		-	
1012001-CCB1		QC		Calibration Blank		-	
C101104-22 A		WC - Alkalinity		SLPO02	No Tag Prefix-A		
C101104-25 A		WC - Alkalinity		SLPO03	No Tag Prefix-A		
C101104-28 A		WC - Alkalinity		SLPO04	No Tag Prefix-A		
C101104-31 A		WC - Alkalinity		SLPO05	No Tag Prefix-A		
C101104-34 A		WC - Alkalinity		SLSW01	No Tag Prefix-A		
C101104-37 A		WC - Alkalinity		SLSW02	No Tag Prefix-A		
C101104-40 A		WC - Alkalinity		SLSW03	No Tag Prefix-A		
C101104-43 A		WC - Alkalinity		SLSW04	No Tag Prefix-A		
C101104-46 A		WC - Alkalinity		SLSW05	No Tag Prefix-A		
C101104-49 A		WC - Alkalinity		SLSWDR3	No Tag Prefix-A		
1012001-CCV2		QC	9072201	Calibration Check		-	
1012001-CCB2		QC		Calibration Blank		-	

ANALYSIS SEQUENCE

1012002

11/30/10

Instrument: Mettler AT

Sequence Date: 11/30/10 00:00

Printed: 12/1/2010 9:58:39AM

Lab Number	Dilut. Factor	Analysis	STD ID	Sample/Std Name	EPA Tag ID	Source Sple	Comments
1011116-SRM1		QC		Reference		-	
1011116-BLK1		QC		Blank		-	
C101104-52 A		WC - Alkalinity		SLSWDR4	No Tag Prefix-A		
1011116-DUP1		QC		Duplicate		C101104-52	
C101104-55 A		WC - Alkalinity		SLSWDR6	No Tag Prefix-A		
C101104-58 A		WC - Alkalinity		SLSWDR7b	No Tag Prefix-A		
C101104-61 A		WC - Alkalinity		SLSWDR7c	No Tag Prefix-A		
C101104-64 A		WC - Alkalinity		SLSWFB	No Tag Prefix-A		
C101104-67 A		WC - Alkalinity		SLSWP06	No Tag Prefix-A		
C101104-70 A		WC - Alkalinity		SLSWP07a	No Tag Prefix-A		
1012002-CCV1		QC	9072201	Calibration Check		-	
1012002-CCB1		QC		Calibration Blank		-	
C101104-73 A		WC - Alkalinity		SLSWP07b	No Tag Prefix-A		
C101104-76 A		WC - Alkalinity		SLSWP08	No Tag Prefix-A		
C101104-79 A		WC - Alkalinity		SLSWP09	No Tag Prefix-A		
C101104-82 A		WC - Alkalinity		SLSWP10	No Tag Prefix-A		
C101104-85 A		WC - Alkalinity		SLSWP11	No Tag Prefix-A		
C101104-88 A		WC - Alkalinity		SLSWP12	No Tag Prefix-A		
C101104-91 A		WC - Alkalinity		SLSWP14	No Tag Prefix-A		
C101104-94 A		WC - Alkalinity		SLSWP15	No Tag Prefix-A		
C101104-97 A		WC - Alkalinity		SLSWPP	No Tag Prefix-A		
1012002-CCV2		QC	9072201	Calibration Check		-	
1012002-CCB2		QC		Calibration Blank		-	

Method: pHCAL Calibration
Start time: 11/30/2010 12:47:50 PM

3/13/2008 2:14:42 PM

11/30/10

Results

No.	Comment / ID	Start time	Rx	Result	Unit	Name
1/2	FISHER (Ref. 25oC) 7.00	11/30/2010 12:47:51 PM	R1	-1.700	mv	--
2/2	FISHER (Ref. 25oC) 4.00	11/30/2010 12:49:29 PM	R1	173.300	mv	--
	Number of segments	1				
	Slope	-58.10 mV/pH				
	Zero point	6.983 pH				
	Calibration temperature	22.0 oC				

Method: Titer Standardization 3/13/2008 2:15:08 PM
Start time: 11/30/2010 1:01:27 PM *lv. 1/30/10*

Results

No.	Comment / ID	Start time	Rx	Result	Unit	Name
1/3	NaOH	11/30/2010 1:01:27 PM	R1	=0.01940	N	Titer
2/3	NaOH	11/30/2010 1:04:48 PM	R1	=0.01990	N	Titer
3/3	NaOH	11/30/2010 1:07:56 PM	R1	=0.01988	N	Titer
-/-	Titer	0.01973	R2	=0.01973	--	Mean Titer

Method: Alkalinity2008 alkalinity
Start time: 11/30/2010 1:12:50 PM

10/16/2009 3:16:38 PM

Results

No.	Comment / ID	Start time	Rx	Result	Unit	Name
1/24	icv/1011115-srm	11/30/2010 1:12:50 PM	R1	=123.480	mg/L	Alkalinity
			R2	=6.17	mL of titrant	ml of titrant
			R3	=8.844	initial pH	initial pH
2/24	icb/1011115-blk	11/30/2010 1:16:00 PM	R1	=0.980	mg/L	Alkalinity
			R2	=0.05	mL of titrant	ml of titrant
			R3	=6.040	initial pH	initial pH
3/24	c101104-01	11/30/2010 1:17:10 PM	R1	=94.700	mg/L	Alkalinity
			R2	=4.74	mL of titrant	ml of titrant
			R3	=5.828	initial pH	initial pH
4/24	1011115-dup	11/30/2010 1:19:53 PM	R1	=95.740	mg/L	Alkalinity
			R2	=4.79	mL of titrant	ml of titrant
			R3	=5.976	initial pH	initial pH
5/24	c101104-04	11/30/2010 1:22:37 PM	R1	=96.180	mg/L	Alkalinity
			R2	=4.81	mL of titrant	ml of titrant
			R3	=5.981	initial pH	initial pH
6/24	c101104-07	11/30/2010 1:25:20 PM	R1	=103.340	mg/L	Alkalinity
			R2	=5.17	mL of titrant	ml of titrant
			R3	=5.923	initial pH	initial pH
7/24	c101104-10	11/30/2010 1:28:10 PM	R1	=104.980	mg/L	Alkalinity
			R2	=5.25	mL of titrant	ml of titrant
			R3	=5.931	initial pH	initial pH
8/24	c101104-13	11/30/2010 1:31:01 PM	R1	=108.660	mg/L	Alkalinity
			R2	=5.43	mL of titrant	ml of titrant
			R3	=5.941	initial pH	initial pH
9/24	c101104-16	11/30/2010 1:33:55 PM	R1	=106.020	mg/L	Alkalinity
			R2	=5.30	mL of titrant	ml of titrant
			R3	=5.964	initial pH	initial pH
10/24	c101104-19	11/30/2010 1:36:49 PM	R1	=133.960	mg/L	Alkalinity
			R2	=6.70	mL of titrant	ml of titrant
			R3	=6.044	initial pH	initial pH

Method: Alkalinity2008 **alkalinity** **10/16/2009 3:16:38 PM**
Start time: 11/30/2010 1:12:50 PM

11/24	ccv1	11/30/2010 1:40:07 PM	R1 =101.900	mg/L	Alkalinity
			R2 =5.10	mL of titrant	ml of titrant
			R3 =9.989	initial pH	initial pH
12/24	ccb1	11/30/2010 1:42:46 PM	R1 =0.580	mg/L	Alkalinity
			R2 =0.03	mL of titrant	ml of titrant
			R3 =5.467	initial pH	initial pH
13/24	c101104-22	11/30/2010 1:43:56 PM	R1 =128.440	mg/L	Alkalinity
			R2 =6.42	mL of titrant	ml of titrant
			R3 =5.820	initial pH	initial pH
14/24	c101104-25	11/30/2010 1:47:11 PM	R1 =122.940	mg/L	Alkalinity
			R2 =6.15	mL of titrant	ml of titrant
			R3 =6.005	initial pH	initial pH
15/24	c101104-28	11/30/2010 1:50:22 PM	R1 =112.100	mg/L	Alkalinity
			R2 =5.60	mL of titrant	ml of titrant
			R3 =5.932	initial pH	initial pH
16/24	c101104-31	11/30/2010 1:53:23 PM	R1 =95.680	mg/L	Alkalinity
			R2 =4.78	mL of titrant	ml of titrant
			R3 =5.921	initial pH	initial pH
17/24	c101104-34	11/30/2010 1:56:09 PM	R1 =131.380	mg/L	Alkalinity
			R2 =6.57	mL of titrant	ml of titrant
			R3 =6.116	initial pH	initial pH
18/24	c101104-37	11/30/2010 1:59:25 PM	R1 =121.860	mg/L	Alkalinity
			R2 =6.09	mL of titrant	ml of titrant
			R3 =6.193	initial pH	initial pH
19/24	c101104-40	11/30/2010 2:02:34 PM	R1 =122.400	mg/L	Alkalinity
			R2 =6.12	mL of titrant	ml of titrant
			R3 =6.257	initial pH	initial pH
20/24	c101104-43	11/30/2010 2:05:44 PM	R1 =122.100	mg/L	Alkalinity
			R2 =6.10	mL of titrant	ml of titrant
			R3 =6.232	initial pH	initial pH
21/24	c101104-46	11/30/2010 2:08:51 PM	R1 =92.820	mg/L	Alkalinity
			R2 =4.64	mL of titrant	ml of titrant
			R3 =5.970	initial pH	initial pH

Method: Alkalinity2008 alkalinity **10/16/2009 3:16:38 PM**
Start time: 11/30/2010 1:12:50 PM

22/24	c101104-49	11/30/2010 2:11:35 PM	R1 =104.960	mg/L	Alkalinity
			R2 =5.25	mL of titrant	ml of titrant
			R3 =5.924	initial pH	initial pH
23/24	ccv2	11/30/2010 2:14:30 PM	R1 =99.780	mg/L	Alkalinity
			R2 =4.99	mL of titrant	ml of titrant
			R3 =10.025	initial pH	initial pH
24/24	ccb2	11/30/2010 2:17:09 PM	R1 =0.940	mg/L	Alkalinity
			R2 =0.05	mL of titrant	ml of titrant
			R3 =5.563	initial pH	initial pH

Method: Alkalinity2008 alkalinity
Start time: 11/30/2010 3:53:53 PM *11/30/10* 10/16/2009 3:16:38 PM

Results

No.	Comment / ID	Start time	Rx	Result	Unit	Name
1/23	icv/1011116-srm	11/30/2010 3:53:53 PM	R1	=124.460	mg/L	Alkalinity
			R2	=6.22	mL of titrant	ml of titrant
			R3	=8.922	initial pH	initial pH
2/23	icb/1011116-blk	11/30/2010 3:57:09 PM	R1	=0.840	mg/L	Alkalinity
			R2	=0.04	mL of titrant	ml of titrant
			R3	=5.930	initial pH	initial pH
3/23	c101104-52	11/30/2010 3:58:20 PM	R1	=103.360	mg/L	Alkalinity
			R2	=5.17	mL of titrant	ml of titrant
			R3	=5.819	initial pH	initial pH
4/23	1011116-dup	11/30/2010 4:01:16 PM	R1	=104.760	mg/L	Alkalinity
			R2	=5.24	mL of titrant	ml of titrant
			R3	=5.903	initial pH	initial pH
5/23	c101104-55	11/30/2010 4:04:09 PM	R1	=130.240	mg/L	Alkalinity
			R2	=6.51	mL of titrant	ml of titrant
			R3	=5.986	initial pH	initial pH
6/23	c101104-58	11/30/2010 4:07:26 PM	R1	=109.960	mg/L	Alkalinity
			R2	=5.50	mL of titrant	ml of titrant
			R3	=6.053	initial pH	initial pH
7/23	c101104-61	11/30/2010 4:10:23 PM	R1	=111.240	mg/L	Alkalinity
			R2	=5.56	mL of titrant	ml of titrant
			R3	=6.050	initial pH	initial pH
8/23	c101104-64	11/30/2010 4:13:21 PM	R1	=1.360	mg/L	Alkalinity
			R2	=0.07	mL of titrant	ml of titrant
			R3	=5.330	initial pH	initial pH
9/23	c101104-67	11/30/2010 4:14:38 PM	R1	=127.040	mg/L	Alkalinity
			R2	=6.35	mL of titrant	ml of titrant
			R3	=5.974	initial pH	initial pH
10/23	c101104-70	11/30/2010 4:17:50 PM	R1	=109.240	mg/L	Alkalinity
			R2	=5.46	mL of titrant	ml of titrant
			R3	=6.117	initial pH	initial pH

Method: Alkalinity2008 **alkalinity** **10/16/2009 3:16:38 PM**
Start time: 11/30/2010 3:53:53 PM

11/23	ccv1	11/30/2010 4:20:48 PM	R1 =100.940	mg/L	Alkalinity
			R2 =5.05	mL of titrant	ml of titrant
12/23	ccb1	11/30/2010 4:23:27 PM	R3 =10.359	initial pH	initial pH
			R1 =1.100	mg/L	Alkalinity
			R2 =0.06	mL of titrant	ml of titrant
13/23	c101104-73	11/30/2010 4:24:39 PM	R3 =5.587	initial pH	initial pH
			R1 =111.380	mg/L	Alkalinity
			R2 =5.57	mL of titrant	ml of titrant
14/23	c101104-76	11/30/2010 4:27:40 PM	R3 =6.003	initial pH	initial pH
			R1 =105.520	mg/L	Alkalinity
			R2 =5.28	mL of titrant	ml of titrant
15/23	c101104-79	11/30/2010 4:30:38 PM	R3 =6.127	initial pH	initial pH
			R1 =101.112	mg/L	Alkalinity
			R2 =5.06	mL of titrant	ml of titrant
16/23	c101104-82	11/30/2010 4:33:32 PM	R3 =6.231	initial pH	initial pH
			R1 =134.420	mg/L	Alkalinity
			R2 =6.72	mL of titrant	ml of titrant
17/23	c101104-85	11/30/2010 4:36:52 PM	R3 =6.071	initial pH	initial pH
			R1 =7.080	mg/L	Alkalinity
			R2 =0.35	mL of titrant	ml of titrant
18/23	c101104-88	11/30/2010 4:38:23 PM	R3 =5.172	initial pH	initial pH
			R1 =99.920	mg/L	Alkalinity
			R2 =5.00	mL of titrant	ml of titrant
19/23	c101104-91	11/30/2010 4:41:13 PM	R3 =6.301	initial pH	initial pH
			R1 =100.180	mg/L	Alkalinity
			R2 =5.01	mL of titrant	ml of titrant
20/23	c101104-94	11/30/2010 4:44:05 PM	R3 =6.303	initial pH	initial pH
			R1 =100.780	mg/L	Alkalinity
			R2 =5.04	mL of titrant	ml of titrant
21/23	c101104-97	11/30/2010 4:46:57 PM	R3 =6.258	initial pH	initial pH
			R1 =105.526	mg/L	Alkalinity
			R2 =5.28	mL of titrant	ml of titrant
			R3 =6.059	initial pH	initial pH

Method:	Alkalinity2008	alkalinity	10/16/2009 3:16:38 PM
Start time:	11/30/2010 3:53:53 PM		

22/23	ccv2	11/30/2010 4:50:08 PM	R1 =99.780	mg/L	Alkalinity
			R2 =4.99	mL of titrant	ml of titrant
			R3 =10.340	initial pH	initial pH
23/23	ccv2	11/30/2010 4:52:47 PM	R1 =0.937	mg/L	Alkalinity
			R2 =0.05	mL of titrant	ml of titrant
			R3 =5.369	initial pH	initial pH

PERKIN ELMER OPTIMA 4300DV ICP-OE

Project(s): Rico-Argentine Water Nov 2010 Date: 11 / 29 / 2010
 Work Order(s): C101104 TDF: DG-220 Analyst: Scott VanOvermelen

Batch Preparation Information

Digest / Prep

TR / Total (Diss)

Matrix

(Water) / Soil / Other

Batch ID

101109 / 101110

Data Storage

Data File: X:/Metals Data Files/

DG220-101109-OED 101129

Standard Information

Calibration Std. # 1 = Reagent Blank Solution

Calibration Std. # 2 = ESAT High, SS-114-039

ICV: SS-114-040 (LIMS ID: 0092802)

Prepped: 09/27/10 By: SWPrepped: 09/27/10 By: SW

CCV: 1:2 of SS-114-039 (LIMS ID: 0092805)

Prepped: 11/29/10 By: SV

CRQL Stock: SS-114-043 (LIMS ID: 0092901)

Prepped: 9-27-2010 By: SW

CRQL Daily (LIMS ID: 0092902)

Prepped: 11/29/10 By: SV

ICSA: SS-114-041 Prepped: 9-27-2010

ICSAB: SS-114-042 Prepped: 9-27-2010

Spike Information

Dissolved Spikes

Sample ID: C101104-03Sample ID: C101104-54Sample vol: 10 mLPSS2007-219: 100 uL

QCS-3 Exp: 5-1-2011 (LIMS ID: 0092734)

SS-114-044: 100 uL

Salt Spike Prepped 9-30-2010 (LIMS ID: 0093002)

Tot. / Tot. Rec. Spikes

Sample ID: -Sample Vol: - mLPSS2007-220: - uL

WW-LFS1 Exp: 5-1-2011 (LIMS ID: 0092735)

PSS2007-221: - uL

WW-LFS2 Exp: 5-1-2011 (LIMS ID: 0092736)

Comments / Maintenance

Replace Nebulizer? Y (N)New pump tubing? Y (N)Replace torch or injector? Y (N)

Analytes Reported:

Ag, (Al), As, (Ba), Be, B, (Ca), Cd, Co, Cr, Cu, (Fe), (K), (Mg)Mn, Mo, (Na), Ni, Pb, Sb, Se, SiO₂, Sr, Ti, Tl, V, (Zn)Sequence ID: 101113Lims Entry (Date / Init): 11/29/10 SV

TLF-06.01		SOP: QAQ-04.00		Eff. Date: 1/17/2007	
ESAT Region 8					
Analytical Data Review Form					
Analyst / Bench Review – Level I					
LIMS: C101104		Project: Rico-Argentine Waters Nov 2010			
TDF: DG-220		Matrix: water		Analysis: Dissolved metals	
Method / Instrument QC Parameters			Analytical Batch / Sample Parameters		
<input type="checkbox"/> Yes <input type="checkbox"/> No <i>N/A</i> Corr. Coef. ≥ 0.995		<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No ICV 90% - 110%		<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No Mth. Blk. (MB) / Prep. BLK (PB) $\leq \pm$ PQL	
<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No ICB $\leq \pm$ PQL		<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No CRI 70% - 130%		<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No Blk. Spike (BS) / LCS (SRM) In Control	
<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No ICSA Spiked Analytes 80% - 120%		<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No ICSA Non-Spiked Analytes $\leq \pm$ PQL		<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No Laboratory Duplicate Analyzed	
<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No ICSAB Spiked Analytes 80% - 120%		<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No CCBs $\leq \pm$ PQL		<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No MS / MSD Analyzed	
<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No CCBs $\leq \pm$ PQL		<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No CCVs 90% - 110%		<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No Serial Dilution Analyzed	
<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No QA QC Samples Internal Stds. 60% - 125%		<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No		<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No <i>N/A</i> Post-Digest Spike Analyzed for Tot. / Tot. Rec.	
<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No		<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No		<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No Samples Internal Stds. 60% - 125%	
Other obvious data quality issues are identified <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No					
Describe any anomaly or deficiency not indicated above in the space provided					
LIMS Electronic Data Transfer					
<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No The instrument data file is uploaded to the X: drive			<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No Instrument data are uploaded into the LIMS		
<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No All samples and QC data are present in LIMS			<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No The analyte list for the sequence is complete		
Analyst: <i>[Signature]</i>			Date: 11/29/10		
Peer Review of Analytical Analysis – Level II					
Method / Instrument QC Parameters			Analytical Batch / Sample Parameters		
<input type="checkbox"/> Yes <input type="checkbox"/> No <i>N/A</i> Corr. Coef. ≥ 0.995		<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No ICV 90% - 110%		<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No Mth. Blk. (MB) / Prep. BLK (PB) $\leq \pm$ PQL	
<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No ICB $\leq \pm$ PQL		<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No CRI 70% - 130%		<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No Blk. Spike (BS) / LCS (SRM) In Control	
<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No ICSA Spiked Analytes 80% - 120%		<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No ICSA Non-Spiked Analytes $\leq \pm$ PQL		<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No Laboratory Duplicate Analyzed	
<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No ICSAB Spiked Analytes 80% - 120%		<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No CCBs $\leq \pm$ PQL		<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No MS / MSD Analyzed	
<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No CCBs $\leq \pm$ PQL		<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No CCVs 90% - 110%		<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No Serial Dilution Analyzed	
<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No QA QC Samples Internal Stds. 60% - 125%		<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No		<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No <i>N/A</i> Post-Digest Spike Analyzed for Tot. / Tot. Rec.	
<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No		<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No		<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No Samples Internal Stds. 60% - 125%	
Other obvious data quality issues are identified <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No					
Describe any anomaly or deficiency not indicated above in the space provided					
LIMS Electronic Data Transfer					
<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No The instrument data file is uploaded to the X: drive			<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No Instrument data are uploaded into the LIMS		
<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No All samples and QC data are present in LIMS			<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No The analyte list for the sequence is complete		
Peer Reviewer: <i>[Signature]</i>			Date: 12/2/2010		

PREPARATION BENCH SHEET

1011109

Matrix: Water

Date Prepared: 11/24/10 11:24 By: SV

TechLaw, Inc. - ESAT Region 8

Prepared using: ADMIN - No Lab Prep Reqd

Printed: 11/24/2010 11:25:36AM

Lab Number	Analysis	EPA Tag ID	Initial (mL)	Final (mL)	Spike1 ID	ul Spike1	Spike2 ID	ul Spike2	Source ID	QC Code	Extraction Comments
C101104-03 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLDRBG	
C101104-06 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLDRBG DUP	
C101104-09 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLDRMZ1a	
C101104-12 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLDRMZ1b	
C101104-15 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLDRMZ1c	
C101104-18 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLDRMZ2	
C101104-21 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLPO01	
C101104-24 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLPO02	
C101104-27 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLPO03	
C101104-30 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLPO04	
C101104-33 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLPO05	
C101104-36 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLSW01	
C101104-39 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLSW02	
C101104-42 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLSW03	
C101104-45 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLSW04	
C101104-48 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLSW05	
C101104-51 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLSWDR3	
C101104-03 A	ICPOE Diss. Metals-2010	lo Tag Prefix-	50	50						SLDRBG	
C101104-06 A	ICPOE Diss. Metals-2010	lo Tag Prefix-	50	50						SLDRBG DUP	
C101104-09 A	ICPOE Diss. Metals-2010	lo Tag Prefix-	50	50						SLDRMZ1a	
C101104-12 A	ICPOE Diss. Metals-2010	lo Tag Prefix-	50	50						SLDRMZ1b	
C101104-15 A	ICPOE Diss. Metals-2010	lo Tag Prefix-	50	50						SLDRMZ1c	
C101104-18 A	ICPOE Diss. Metals-2010	lo Tag Prefix-	50	50						SLDRMZ2	

PREPARATION BENCH SHEET

1011109

TechLaw, Inc. - ESAT Region 8

Printed: 11/24/2010 11:25:36AM

Matrix: Water

Date Prepared: 11/24/10 11:24 By: SV

Prepared using: ADMIN - No Lab Prep Req'd

Lab Number	Analysis	EPA Tag ID	Initial (mL)	Final (mL)	Spike 1 ID	ul Spike1	Spike2 ID	ul Spike2	Source ID	QC Code	Extraction Comments
C101104-21 A	ICPOE Diss. Metals-2010	IC Tag Prefix-	50	50						SLPO01	
C101104-24 A	ICPOE Diss. Metals-2010	IC Tag Prefix-	50	50						SLPO02	
C101104-27 A	ICPOE Diss. Metals-2010	IC Tag Prefix-	50	50						SLPO03	
C101104-30 A	ICPOE Diss. Metals-2010	IC Tag Prefix-	50	50						SLPO04	
C101104-33 A	ICPOE Diss. Metals-2010	IC Tag Prefix-	50	50						SLPO05	
C101104-36 A	ICPOE Diss. Metals-2010	IC Tag Prefix-	50	50						SLSW01	
C101104-39 A	ICPOE Diss. Metals-2010	IC Tag Prefix-	50	50						SLSW02	
C101104-42 A	ICPOE Diss. Metals-2010	IC Tag Prefix-	50	50						SLSW03	
C101104-45 A	ICPOE Diss. Metals-2010	IC Tag Prefix-	50	50						SLSW04	
C101104-48 A	ICPOE Diss. Metals-2010	IC Tag Prefix-	50	50						SLSW05	
C101104-51 A	ICPOE Diss. Metals-2010	IC Tag Prefix-	50	50						SLSWDR3	
1011109-BLK1	QC		50	50						Blank	
1011109-BS1	QC		10	10	0100114	100				LCS	
1011109-DUP1	QC		50	50					C101104-03	Duplicate	
1011109-MS1	QC		10	10	0100114	100			C101104-03	Matrix Spike	
1011109-MSD1	QC		10	10	0100114	100			C101104-03	Matrix Spike Dup	

PREPARATION BENCH SHEET

1011110

TechLaw, Inc. - ESAT Region 8

Matrix: Water

Date Prepared: 11/24/10 11:25 By: SV

Prepared using: ADMIN - No Lab Prep Req'd

Lab Number	Analysis	EPA Tag ID	Initial (mL)	Final (mL)	Spike1 ID	ul Spike1	Spike2 ID	ul Spike2	Source ID	QC Code	Extraction Comments
C101104-54 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLSWDR4	
C101104-57 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLSWDR6	
C101104-60 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLSWDR7b	
C101104-63 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLSWDR7c	
C101104-66 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLSWFB	
C101104-69 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLSWP06	
C101104-72 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLSWP07a	
C101104-75 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLSWP07b	
C101104-78 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLSWP08	
C101104-81 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLSWP09	
C101104-84 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLSWP10	
C101104-87 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLSWP11	
C101104-90 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLSWP12	
C101104-93 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLSWP14	
C101104-96 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLSWP15	
C101104-99 A	DM-Hardness - Calculated	lo Tag Prefix-	50	50						SLSWPP	
C101104-54 A	ICPOE Diss. Metals-2010	lo Tag Prefix-	50	50						SLSWDR4	
C101104-57 A	ICPOE Diss. Metals-2010	lo Tag Prefix-	50	50						SLSWDR6	
C101104-60 A	ICPOE Diss. Metals-2010	lo Tag Prefix-	50	50						SLSWDR7b	
C101104-63 A	ICPOE Diss. Metals-2010	lo Tag Prefix-	50	50						SLSWDR7c	
C101104-66 A	ICPOE Diss. Metals-2010	lo Tag Prefix-	50	50						SLSWFB	
C101104-69 A	ICPOE Diss. Metals-2010	lo Tag Prefix-	50	50						SLSWP06	
C101104-72 A	ICPOE Diss. Metals-2010	lo Tag Prefix-	50	50						SLSWP07a	

PREPARATION BENCH SHEET

1011110

Matrix: Water

TechLaw, Inc. - ESAT Region 8

Printed: 11/24/2010 11:26:55AM

Date Prepared: 11/24/10 11:25 By: SV

Prepared using: ADMIN - No Lab Prep Req'd

Lab Number	Analysis	EPA Tag ID	Initial (mL)	Final (mL)	Spike1 ID	ul Spike1	Spike2 ID	ul Spike2	Source ID	QC Code	Extraction Comments
C101104-75 A	ICPOE Diss. Metals-2010	IC Tag Prefix-	50	50						SLSWP07b	
C101104-78 A	ICPOE Diss. Metals-2010	IC Tag Prefix-	50	50						SLSWP08	
C101104-81 A	ICPOE Diss. Metals-2010	IC Tag Prefix-	50	50						SLSWP09	
C101104-84 A	ICPOE Diss. Metals-2010	IC Tag Prefix-	50	50						SLSWP10	
C101104-87 A	ICPOE Diss. Metals-2010	IC Tag Prefix-	50	50						SLSWP11	
C101104-90 A	ICPOE Diss. Metals-2010	IC Tag Prefix-	50	50						SLSWP12	
C101104-93 A	ICPOE Diss. Metals-2010	IC Tag Prefix-	50	50						SLSWP14	
C101104-96 A	ICPOE Diss. Metals-2010	IC Tag Prefix-	50	50						SLSWP15	
C101104-99 A	ICPOE Diss. Metals-2010	IC Tag Prefix-	50	50						SLSWPP	
1011110-BLK1	QC		50	50						Blank	
1011110-BS1	QC		10	10	0100114	100				LCS	
1011110-DUP1	QC		50	50					C101104-54	Duplicate	
1011110-MS1	QC		10	10	0100114	100			C101104-54	Matrix Spike	
1011110-MSD1	QC		10	10	0100114	100			C101104-54	Matrix Spike Dup	

ANALYSIS SEQUENCE

1011113

11/29/10

Instrument: ICPOE - PE Optima

Sequence Date: 11/29/10 00:00

Printed: 11/29/2010 11:11:35AM

Lab Number	Dilut. Factor	Analysis	STD ID	Sample/Std Name	EPA Tag ID	Source Sple	Comments
1011113-ICV1		QC	0092802	Initial Cal Check		-	
1011113-ICB1		QC		Initial Cal Blank		-	
1011113-CRL1		QC	0092902	MRL Check		-	
1011113-IFA1		QC	0092803	Interference Check A		-	
1011113-IFB1		QC	0092804	Interference Check B		-	
1011113-CCV1		QC	0092805	Calibration Check		-	
1011113-CCB1		QC		Calibration Blank		-	
1011109-BLK1		QC		Blank		-	
1011109-BS1		QC		LCS		-	
C101104-03 A		DM-Hardness - Calculated		SLDRBG	No Tag Prefix-C		
C101104-03 A		ICPOE Diss. Metals-2010		SLDRBG	No Tag Prefix-C		
1011109-DUP1		QC		Duplicate		C101104-03	
1011113-SRD1		QC		Serial Dilution		C101104-03	
1011109-MS1		QC		Matrix Spike		C101104-03	
1011109-MSD1		QC		Matrix Spike Dup		C101104-03	
C101104-06 A		DM-Hardness - Calculated		SLDRBG DUP	No Tag Prefix-C		
C101104-06 A		ICPOE Diss. Metals-2010		SLDRBG DUP	No Tag Prefix-C		
C101104-09 A		DM-Hardness - Calculated		SLDRMZ1a	No Tag Prefix-C		
C101104-09 A		ICPOE Diss. Metals-2010		SLDRMZ1a	No Tag Prefix-C		
1011113-CCV2		QC	0092805	Calibration Check		-	
1011113-CCB2		QC		Calibration Blank		-	
C101104-12 A		DM-Hardness - Calculated		SLDRMZ1b	No Tag Prefix-C		
C101104-12 A		ICPOE Diss. Metals-2010		SLDRMZ1b	No Tag Prefix-C		
C101104-15 A		DM-Hardness - Calculated		SLDRMZ1c	No Tag Prefix-C		
C101104-15 A		ICPOE Diss. Metals-2010		SLDRMZ1c	No Tag Prefix-C		
C101104-18 A		DM-Hardness - Calculated		SLDRMZ2	No Tag Prefix-C		
C101104-18 A		ICPOE Diss. Metals-2010		SLDRMZ2	No Tag Prefix-C		
C101104-21 A		DM-Hardness - Calculated		SLPO01	No Tag Prefix-C		
C101104-21 A		ICPOE Diss. Metals-2010		SLPO01	No Tag Prefix-C		
C101104-24 A		DM-Hardness - Calculated		SLPO02	No Tag Prefix-C		
C101104-24 A		ICPOE Diss. Metals-2010		SLPO02	No Tag Prefix-C		
C101104-27 A		DM-Hardness - Calculated		SLPO03	No Tag Prefix-C		
C101104-27 A		ICPOE Diss. Metals-2010		SLPO03	No Tag Prefix-C		
C101104-30 A		DM-Hardness - Calculated		SLPO04	No Tag Prefix-C		
C101104-30 A		ICPOE Diss. Metals-2010		SLPO04	No Tag Prefix-C		
C101104-33 A		DM-Hardness - Calculated		SLPO05	No Tag Prefix-C		
C101104-33 A		ICPOE Diss. Metals-2010		SLPO05	No Tag Prefix-C		
C101104-36 A		DM-Hardness - Calculated		SLSW01	No Tag Prefix-C		
C101104-36 A		ICPOE Diss. Metals-2010		SLSW01	No Tag Prefix-C		
1011113-CCV3		QC	0092805	Calibration Check		-	

ANALYSIS SEQUENCE

1011113

Instrument: ICPOE - PE Optima

Sequence Date: 11/29/10 00:00

Printed: 11/29/2010 11:11:35AM

Lab Number	Dilut. Factor	Analysis	STD ID	Sample/Std Name	EPA Tag ID	Source Sple	Comments
1011113-CCB3		QC		Calibration Blank		-	
C101104-39 A		DM-Hardness - Calculated		SLSW02	No Tag Prefix-C		
C101104-39 A		ICPOE Diss. Metals-2010		SLSW02	No Tag Prefix-C		
C101104-42 A		DM-Hardness - Calculated		SLSW03	No Tag Prefix-C		
C101104-42 A		ICPOE Diss. Metals-2010		SLSW03	No Tag Prefix-C		
C101104-45 A		DM-Hardness - Calculated		SLSW04	No Tag Prefix-C		
C101104-45 A		ICPOE Diss. Metals-2010		SLSW04	No Tag Prefix-C		
C101104-48 A		DM-Hardness - Calculated		SLSW05	No Tag Prefix-C		
C101104-48 A		ICPOE Diss. Metals-2010		SLSW05	No Tag Prefix-C		
C101104-51 A		DM-Hardness - Calculated		SLSWDR3	No Tag Prefix-C		
C101104-51 A		ICPOE Diss. Metals-2010		SLSWDR3	No Tag Prefix-C		
1011110-BLK1		QC		Blank		-	
1011110-BS1		QC		LCS		-	
C101104-54 A		DM-Hardness - Calculated		SLSWDR4	No Tag Prefix-C		
C101104-54 A		ICPOE Diss. Metals-2010		SLSWDR4	No Tag Prefix-C		
1011110-DUP1		QC		Duplicate		C101104-54	
1011113-CCV4		QC	0092805	Calibration Check		-	
1011113-CCB4		QC		Calibration Blank		-	
1011113-SRD2		QC		Serial Dilution		C101104-54	
1011110-MS1		QC		Matrix Spike		C101104-54	
1011110-MSD1		QC		Matrix Spike Dup		C101104-54	
C101104-57 A		DM-Hardness - Calculated		SLSWDR6	No Tag Prefix-C		
C101104-57 A		ICPOE Diss. Metals-2010		SLSWDR6	No Tag Prefix-C		
C101104-60 A		DM-Hardness - Calculated		SLSWDR7b	No Tag Prefix-C		
C101104-60 A		ICPOE Diss. Metals-2010		SLSWDR7b	No Tag Prefix-C		
C101104-63 A		DM-Hardness - Calculated		SLSWDR7c	No Tag Prefix-C		
C101104-63 A		ICPOE Diss. Metals-2010		SLSWDR7c	No Tag Prefix-C		
C101104-66 A		DM-Hardness - Calculated		SLSWFB	No Tag Prefix-C		
C101104-66 A		ICPOE Diss. Metals-2010		SLSWFB	No Tag Prefix-C		
C101104-69 A		DM-Hardness - Calculated		SLSWP06	No Tag Prefix-C		
C101104-69 A		ICPOE Diss. Metals-2010		SLSWP06	No Tag Prefix-C		
C101104-72 A		DM-Hardness - Calculated		SLSWP07a	No Tag Prefix-C		
C101104-72 A		ICPOE Diss. Metals-2010		SLSWP07a	No Tag Prefix-C		
1011113-CCV5		QC	0092805	Calibration Check		-	
1011113-CCB5		QC		Calibration Blank		-	
C101104-75 A		DM-Hardness - Calculated		SLSWP07b	No Tag Prefix-C		
C101104-75 A		ICPOE Diss. Metals-2010		SLSWP07b	No Tag Prefix-C		
C101104-78 A		DM-Hardness - Calculated		SLSWP08	No Tag Prefix-C		
C101104-78 A		ICPOE Diss. Metals-2010		SLSWP08	No Tag Prefix-C		
C101104-81 A		DM-Hardness - Calculated		SLSWP09	No Tag Prefix-C		

ANALYSIS SEQUENCE

1011113

Instrument: ICPOE - PE Optima

Sequence Date: 11/29/10 00:00

Printed: 11/29/2010 11:11:35AM

Lab Number	Dilut. Factor	Analysis	STD ID	Sample/Std Name	EPA Tag ID	Source Sple	Comments
C101104-81 A		ICPOE Diss. Metals-2010		SLSWP09	No Tag Prefix-C		
C101104-84 A		DM-Hardness - Calculated		SLSWP10	No Tag Prefix-C		
C101104-84 A		ICPOE Diss. Metals-2010		SLSWP10	No Tag Prefix-C		
C101104-87 A		DM-Hardness - Calculated		SLSWP11	No Tag Prefix-C		
C101104-87 A		ICPOE Diss. Metals-2010		SLSWP11	No Tag Prefix-C		
C101104-90 A		DM-Hardness - Calculated		SLSWP12	No Tag Prefix-C		
C101104-90 A		ICPOE Diss. Metals-2010		SLSWP12	No Tag Prefix-C		
C101104-93 A		DM-Hardness - Calculated		SLSWP14	No Tag Prefix-C		
C101104-93 A		ICPOE Diss. Metals-2010		SLSWP14	No Tag Prefix-C		
C101104-96 A		DM-Hardness - Calculated		SLSWP15	No Tag Prefix-C		
C101104-96 A		ICPOE Diss. Metals-2010		SLSWP15	No Tag Prefix-C		
C101104-99 A		DM-Hardness - Calculated		SLSWPP	No Tag Prefix-C		
C101104-99 A		ICPOE Diss. Metals-2010		SLSWPP	No Tag Prefix-C		
1011113-CCV6		QC	0092805	Calibration Check		-	
1011113-CCB6		QC		Calibration Blank		-	

Analytical Sequence

Method : ESAT_2009_1.1

11/11/29/10

Seq.	Loc.		Sample ID
1	1		Cal Blank
2	9		High Std
3	10		SEQ-ICV
4	1		SEQ-ICB
5	11		SEQ-CRL
6	12		SEQ-IFA
7	13		SEQ-IFB
8	3		SEQ-CCV 1
9	1		SEQ-CCB 1
10	26		1011109-BLK1
11	27		1011109-BS1
12	28		C101104-03
13	29		1011109-DUP1
14	30		SEQ-SRD1 @5X
15	31		1011109-MS1
16	32		1011109-MSD1
17	33		C101104-06
18	34		C101104-09
19	35		Blank
20	3		SEQ-CCV 2
21	1		SEQ-CCB 2
22	36		C101104-12
23	37		C101104-15
24	38		C101104-18
25	39		C101104-21
26	40		C101104-24
27	41		C101104-27
28	42		C101104-30
29	43		C101104-33
30	44		C101104-36
31	45		Blank
32	3		SEQ-CCV 3
33	1		SEQ-CCB 3
34	46		C101104-39
35	47		C101104-42
36	48		C101104-45
37	49		C101104-48
38	50		C101104-51
39	51		1011110-BLK1
40	52		1011110-BS1
41	53		C101104-54
42	54		1011110-DUP1
43	55		Blank
44	3		SEQ-CCV 4
45	1		SEQ-CCB 4
46	56		SEQ-SRD1 @5X
47	57		1011110-MS1
48	58		1011110-MSD1
49	59		C101104-57
50	60		C101104-60
51	61		C101104-63
52	62		C101104-66
53	63		C101104-69
54	64		C101104-72
55	65		Blank
56	3		SEQ-CCV 5

Analytical Sequence

Method : ESAT_2009_1.1

Seq.	Loc.		Sample ID
57	1	QC	SEQ-CCB 5
58	66		C101104-75
59	67		C101104-78
60	68		C101104-81
61	69		C101104-84
62	70		C101104-87
63	71		C101104-90
64	72		C101104-93
65	73		C101104-96
66	74		C101104-99
67	75		Blank
68	3	QC	SEQ-CCV 6
69	1	QC	SEQ-CCB 6

Sample Information Detail Report
Document Name: DG220_1011109_OED_101129

File Description

DG-220 Rico-Argentina Waters 2010

Parameters Common to All Samples

Batch ID	1011109/1011110
Analyst Name	Walker
Volume Units	mL
Weight Units	g

for 11/29/10

Parameters That Vary By Sample

Sample No	A/S Location	Sample ID	Remarks
1	26	1011109-BLK1	
2	27	1011109-BS1	
3	28	C101104-03	
4	29	1011109-DUP1	
5	30	SEQ-SRD1 @5X	
6	31	1011109-MS1	
7	32	1011109-MSD1	
8	33	C101104-06	
9	34	C101104-09	
10	35	Blank	
11	36	C101104-12	
12	37	C101104-15	
13	38	C101104-18	
14	39	C101104-21	
15	40	C101104-24	
16	41	C101104-27	
17	42	C101104-30	
18	43	C101104-33	
19	44	C101104-36	
20	45	Blank	
21	46	C101104-39	
22	47	C101104-42	
23	48	C101104-45	
24	49	C101104-48	
25	50	C101104-51	
26	51	1011110-BLK1	
27	52	1011110-BS1	
28	53	C101104-54	
29	54	1011110-DUP1	
30	55	Blank	
31	56	SEQ-SRD1 @5X	
32	57	1011110-MS1	
33	58	1011110-MSD1	
34	59	C101104-57	
35	60	C101104-60	
36	61	C101104-63	
37	62	C101104-66	
38	63	C101104-69	
39	64	C101104-72	
40	65	Blank	
41	66	C101104-75	
42	67	C101104-78	
43	68	C101104-81	
44	69	C101104-84	
45	70	C101104-87	
46	71	C101104-90	
47	72	C101104-93	
48	73	C101104-96	
49	74	C101104-99	
50	75	Blank	

Sample No	Aliquot Volume	Diluted To Vol.	Matrix Check Sample
1			
2			
3			

Recovery 3 of 1

Sample Information Detail Report
Document Name: DG220_1011109_OED_101129

4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50

2

10

Duplicate of 3
5X Dilution of 3
Recovery 3 of 3
Recovery 3 of 3

Recovery 3 of 26

Duplicate of 28

2

10

5X Dilution of 28
Recovery 3 of 28
Recovery 3 of 28

=====

Analysis Begun

Start Time: 11/29/2010 7:37:57 AM

Plasma On Time: 11/29/2010 6:29:07 AM

Logged In Analyst: esat

Technique: ICP Continuous

Spectrometer Model: Optima 4300 DV, S/N 077N3082602 Autosampler Model: AS-93plus

Sample Information File: C:\pe\Administrator\Sample Information\2010\DG220_1011109_OED_101129.sif

Batch ID: 1011109/1011110

Results Data Set: DG220_1011109_101129

Results Library: C:\pe\Administrator\Results\Results.mdb

Sequence No.: 1

Autosampler Location: 1

Sample ID: Cal Blank

Date Collected: 11/29/2010 7:37:58 AM

Analyst:

Data Type: Original

Initial Sample Wt:

Initial Sample Vol:

Dilution:

Sample Prep Vol:

Nebulizer Parameters: Cal Blank

Analyte	Back Pressure	Flow
All	115.0 kPa	0.80 L/min

Mean Data: Cal Blank

Analyte	Mean Corrected Intensity	Std.Dev.	RSD	Conc.	Calib Units
Sc Axial	2450827.2	10476.22	0.43%	100.00	%
Sc Radial	299727.1	762.79	0.25%	100.0	%
Ag 328.068†	-225.6	11.70	5.19%	[0.00]	ug/L
Al 396.153†	32.7	11.34	34.65%	[0.00]	ug/L
As 193.696†	-7.5	1.43	18.91%	[0.00]	ug/L
Ba 233.527†	-18.2	1.24	6.81%	[0.00]	ug/L
Be 313.107†	-120.3	8.11	6.74%	[0.00]	ug/L
B 249.677†	-303.5	5.11	1.68%	[0.00]	ug/L
Ca 317.933†	1.3	1.02	77.29%	[0.00]	ug/L
Cd 214.440†	-12.5	1.04	8.28%	[0.00]	ug/L
Co 228.616†	-42.0	3.04	7.24%	[0.00]	ug/L
Cr 267.716†	6.4	0.44	6.78%	[0.00]	ug/L
Cu 324.752†	550.5	11.19	2.03%	[0.00]	ug/L
Fe 238.204†	11.9	0.94	7.87%	[0.00]	ug/L
K 766.490†	350.5	35.02	9.99%	[0.00]	ug/L
Mg 285.213†	20.6	2.79	13.57%	[0.00]	ug/L
Mn 257.610†	35.4	6.86	19.40%	[0.00]	ug/L
Mo 202.031†	-5.1	0.58	11.56%	[0.00]	ug/L
Na 589.592†	209.9	20.06	9.56%	[0.00]	ug/L
Ni 231.604†	28.6	0.27	0.96%	[0.00]	ug/L
Pb 220.353†	-1.0	2.80	283.75%	[0.00]	ug/L
Sb 206.836†	12.9	1.21	9.36%	[0.00]	ug/L
Se 196.026†	-0.1	4.02	>999.9%	[0.00]	ug/L
SiO2 251.603†	162.8	8.68	5.33%	[0.00]	ug/L
Sr 421.552†	-28.3	15.81	55.91%	[0.00]	ug/L
Ti 334.940†	14.6	4.42	30.20%	[0.00]	ug/L
Tl 190.801†	-5.7	1.21	21.06%	[0.00]	ug/L
V 290.880†	2315.4	15.45	0.67%	[0.00]	ug/L
Zn 206.200†	-0.9	1.26	133.84%	[0.00]	ug/L

Sequence No.: 2
 Sample ID: High Std
 Analyst:
 Initial Sample Wt:
 Dilution:

Autosampler Location: 9
 Date Collected: 11/29/2010 7:41:01 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: High Std

Analyte	Back Pressure	Flow
All	115.0 kPa	0.80 L/min

Mean Data: High Std

Analyte	Mean Corrected Intensity	Std.Dev.	RSD	Conc. Units	Calib
Sc Axial	2430883.0	8687.95	0.36%	99.19	%
Sc Radial	303406.5	1879.93	0.62%	101.2	%
Ag 328.068†	30727.1	153.28	0.50%	[500]	ug/L
Al 396.153†	116443.3	707.41	0.61%	[25000]	ug/L
As 193.696†	495.1	3.61	0.73%	[5000]	ug/L
Ba 233.527†	11965.4	85.65	0.72%	[1000]	ug/L
Be 313.107†	274045.8	1355.44	0.49%	[1000]	ug/L
B 249.677†	65122.9	672.40	1.03%	[10000]	ug/L
Ca 317.933†	72046.9	440.00	0.61%	[25000]	ug/L
Cd 214.440†	3275.6	11.96	0.37%	[1000]	ug/L
Co 228.616†	4253.8	16.70	0.39%	[1000]	ug/L
Cr 267.716†	31550.8	193.01	0.61%	[5000]	ug/L
Cu 324.752†	192036.4	138.11	0.07%	[2000]	ug/L
Fe 238.204†	3502.0	36.77	1.05%	[25000]	ug/L
K 766.490†	87938.3	306.40	0.35%	[50000]	ug/L
Mg 285.213†	160474.0	918.42	0.57%	[25000]	ug/L
Mn 257.610†	273982.2	2092.45	0.76%	[2000]	ug/L
Mo 202.031†	492.3	1.01	0.21%	[1000]	ug/L
Na 589.592†	136069.1	490.41	0.36%	[25000]	ug/L
Ni 231.604†	9025.7	83.20	0.92%	[5000]	ug/L
Pb 220.353†	2214.0	4.07	0.18%	[5000]	ug/L
Sb 206.836†	1297.8	3.17	0.24%	[5000]	ug/L
Se 196.026†	424.2	4.08	0.96%	[5000]	ug/L
SiO2 251.603†	66283.2	296.56	0.45%	[20000]	ug/L
Sr 421.552†	1526631.8	3304.58	0.22%	[1000]	ug/L
Ti 334.940†	211398.2	1311.36	0.62%	[1000]	ug/L
Tl 190.801†	1224.2	3.13	0.26%	[5000]	ug/L
V 290.880†	56636.9	419.27	0.74%	[2000]	ug/L
Zn 206.200†	5291.5	69.96	1.32%	[5000]	ug/L

Calibration Summary

Analyte	Stds.	Equation	Intercept	Slope	Curvature	Corr. Coef.	Reslope
Ag 328.068	1	Lin Thru 0	0.0	61.45	0.00000	1.000000	
Al 396.153	1	Lin, Calc Int	0.0	4.658	0.00000	1.000000	
As 193.696	1	Lin Thru 0	0.0	0.0990	0.00000	1.000000	
Ba 233.527	1	Lin Thru 0	0.0	11.97	0.00000	1.000000	
Be 313.107	1	Lin Thru 0	0.0	274.0	0.00000	1.000000	
B 249.677	1	Lin Thru 0	0.0	6.512	0.00000	1.000000	
Ca 317.933	1	Lin Thru 0	0.0	2.882	0.00000	1.000000	
Cd 214.440	1	Lin Thru 0	0.0	3.276	0.00000	1.000000	
Co 228.616	1	Lin Thru 0	0.0	4.254	0.00000	1.000000	
Cr 267.716	1	Lin Thru 0	0.0	6.310	0.00000	1.000000	
Cu 324.752	1	Lin Thru 0	0.0	96.02	0.00000	1.000000	
Fe 238.204	1	Lin, Calc Int	0.0	0.1401	0.00000	1.000000	
K 766.490	1	Lin Thru 0	0.0	1.759	0.00000	1.000000	
Mg 285.213	1	Lin, Calc Int	0.0	6.419	0.00000	1.000000	
Mn 257.610	1	Lin Thru 0	0.0	137.0	0.00000	1.000000	
Mo 202.031	1	Lin Thru 0	0.0	0.4923	0.00000	1.000000	
Na 589.592	1	Lin, Calc Int	0.0	5.443	0.00000	1.000000	
Ni 231.604	1	Lin Thru 0	0.0	1.805	0.00000	1.000000	
Pb 220.353	1	Lin Thru 0	0.0	0.4428	0.00000	1.000000	
Sb 206.836	1	Lin Thru 0	0.0	0.2596	0.00000	1.000000	

Se 196.026	1	Lin Thru 0	0.0	0.0848	0.00000	1.000000
SiO2 251.603	1	Lin, Calc Int	0.0	3.314	0.00000	1.000000
Sr 421.552	1	Lin, Calc Int	0.0	1527	0.00000	1.000000
Ti 334.940	1	Lin Thru 0	0.0	211.4	0.00000	1.000000
Tl 190.801	1	Lin Thru 0	0.0	0.2448	0.00000	1.000000
V 290.880	1	Lin Thru 0	0.0	28.32	0.00000	1.000000
Zn 206.200	1	Lin Thru 0	0.0	1.058	0.00000	1.000000

Sequence No.: 3
 Sample ID: SEQ-ICV
 Analyst:
 Initial Sample Wt:
 Dilution:

Autosampler Location: 10
 Date Collected: 11/29/2010 7:44:09 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: SEQ-ICV

Analyte	Back Pressure	Flow
All	117.0 kPa	0.80 L/min

Mean Data: SEQ-ICV

Analyte	Mean Corrected Intensity	Calib Conc. Units	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2509142.7	102.4 %	0.42			0.41%
Sc Radial	306225.1	102.2 %	0.69			0.67%
Ag 328.068†	15587.7	257.2 ug/L	0.37	257.2 ug/L	0.37	0.14%
QC value within limits for Ag 328.068		Recovery = 102.88%				
Al 396.153†	4779.7	1004 ug/L	12.5	1004 ug/L	12.5	1.25%
QC value within limits for Al 396.153		Recovery = 100.42%				
As 193.696†	194.3	1981 ug/L	16.6	1981 ug/L	16.6	0.84%
QC value within limits for As 193.696		Recovery = 99.07%				
Ba 233.527†	12049.7	1005 ug/L	1.6	1005 ug/L	1.6	0.16%
QC value within limits for Ba 233.527		Recovery = 100.51%				
Be 313.107†	271979.8	992.3 ug/L	1.65	992.3 ug/L	1.65	0.17%
QC value within limits for Be 313.107		Recovery = 99.23%				
B 249.677†	6934.5	1065 ug/L	4.7	1065 ug/L	4.7	0.44%
QC value within limits for B 249.677		Recovery = 106.48%				
Ca 317.933†	2906.9	955.3 ug/L	4.84	955.3 ug/L	4.84	0.51%
QC value within limits for Ca 317.933		Recovery = 95.53%				
Cd 214.440†	3343.9	1021 ug/L	3.6	1021 ug/L	3.6	0.35%
QC value within limits for Cd 214.440		Recovery = 102.08%				
Co 228.616†	4328.2	1018 ug/L	2.9	1018 ug/L	2.9	0.28%
QC value within limits for Co 228.616		Recovery = 101.84%				
Cr 267.716†	6306.0	1000 ug/L	4.1	1000 ug/L	4.1	0.41%
QC value within limits for Cr 267.716		Recovery = 100.04%				
Cu 324.752†	95071.9	991.9 ug/L	3.11	991.9 ug/L	3.11	0.31%
QC value within limits for Cu 324.752		Recovery = 99.19%				
Fe 238.204†	136.7	962.1 ug/L	10.38	962.1 ug/L	10.38	1.08%
QC value within limits for Fe 238.204		Recovery = 96.21%				
K 766.490†	9000.2	4874 ug/L	42.8	4874 ug/L	42.8	0.88%
QC value within limits for K 766.490		Recovery = 97.49%				
Mg 285.213†	6601.3	1019 ug/L	7.3	1019 ug/L	7.3	0.71%
QC value within limits for Mg 285.213		Recovery = 101.91%				
Mn 257.610†	139572.9	1018 ug/L	1.2	1018 ug/L	1.2	0.12%
QC value within limits for Mn 257.610		Recovery = 101.80%				
Mo 202.031†	496.2	1008 ug/L	4.3	1008 ug/L	4.3	0.43%
QC value within limits for Mo 202.031		Recovery = 100.79%				
Na 589.592†	5704.9	965.7 ug/L	15.19	965.7 ug/L	15.19	1.57%
QC value within limits for Na 589.592		Recovery = 96.57%				
Ni 231.604†	1833.3	1018 ug/L	7.6	1018 ug/L	7.6	0.75%
QC value within limits for Ni 231.604		Recovery = 101.75%				
Pb 220.353†	898.3	2022 ug/L	8.6	2022 ug/L	8.6	0.43%
QC value within limits for Pb 220.353		Recovery = 101.08%				
Sb 206.836†	524.6	2015 ug/L	7.8	2015 ug/L	7.8	0.39%
QC value within limits for Sb 206.836		Recovery = 100.76%				
Se 196.026†	85.7	1012 ug/L	30.0	1012 ug/L	30.0	2.96%
QC value within limits for Se 196.026		Recovery = 101.16%				
SiO2 251.603†	16499.9	4923 ug/L	9.5	4923 ug/L	9.5	0.19%
QC value within limits for SiO2 251.603		Recovery = 98.46%				
Sr 421.552†	1561146.0	1023 ug/L	3.3	1023 ug/L	3.3	0.32%
QC value within limits for Sr 421.552		Recovery = 102.26%				
Ti 334.940†	212146.3	1004 ug/L	0.8	1004 ug/L	0.8	0.08%
QC value within limits for Ti 334.940		Recovery = 100.35%				
Tl 190.801†	1237.1	5058 ug/L	21.4	5058 ug/L	21.4	0.42%
QC value within limits for Tl 190.801		Recovery = 101.17%				
V 290.880†	28033.6	991.5 ug/L	3.35	991.5 ug/L	3.35	0.34%
QC value within limits for V 290.880		Recovery = 99.15%				

Zn 206.200†	1051.7	986.0 ug/L	8.92	986.0 ug/L	8.92	0.90%
-------------	--------	------------	------	------------	------	-------

QC value within limits for Zn 206.200 Recovery = 98.60%

All analyte(s) passed QC.

Sequence No.: 4
 Sample ID: SEQ-ICB
 Analyst:
 Initial Sample Wt:
 Dilution:

Autosampler Location: 1
 Date Collected: 11/29/2010 7:47:17 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: SEQ-ICB

Analyte Back Pressure Flow
 All 116.0 kPa 0.80 L/min

Mean Data: SEQ-ICB

Analyte	Mean Corrected Intensity	Calib Conc. Units	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2511539.5	102.5 %	0.47			0.45%
Sc Radial	310523.7	103.6 %	1.24			1.20%
Ag 328.068†	8.2	0.1405 ug/L	0.20158	0.1405 ug/L	0.20158	143.45%
QC value within limits for Ag 328.068		Recovery =	Not calculated			
Al 396.153†	-1.6	-0.7355 ug/L	2.63084	-0.7355 ug/L	2.63084	357.70%
QC value within limits for Al 396.153		Recovery =	Not calculated			
As 193.696†	0.3	2.649 ug/L	9.7738	2.649 ug/L	9.7738	368.96%
QC value within limits for As 193.696		Recovery =	Not calculated			
Ba 233.527†	-0.1	-0.0187 ug/L	0.12057	-0.0187 ug/L	0.12057	646.14%
QC value within limits for Ba 233.527		Recovery =	Not calculated			
Be 313.107†	12.7	0.0456 ug/L	0.01965	0.0456 ug/L	0.01965	43.10%
QC value within limits for Be 313.107		Recovery =	Not calculated			
B 249.677†	268.1	41.17 ug/L	2.578	41.17 ug/L	2.578	6.26%
QC value within limits for B 249.677		Recovery =	Not calculated			
Ca 317.933†	1.8	0.5303 ug/L	1.72811	0.5303 ug/L	1.72811	325.87%
QC value within limits for Ca 317.933		Recovery =	Not calculated			
Cd 214.440†	0.6	0.1972 ug/L	0.37193	0.1972 ug/L	0.37193	188.61%
QC value within limits for Cd 214.440		Recovery =	Not calculated			
Co 228.616†	0.3	0.0684 ug/L	0.38559	0.0684 ug/L	0.38559	563.31%
QC value within limits for Co 228.616		Recovery =	Not calculated			
Cr 267.716†	-0.6	-0.0953 ug/L	0.24033	-0.0953 ug/L	0.24033	252.22%
QC value within limits for Cr 267.716		Recovery =	Not calculated			
Cu 324.752†	12.6	0.1338 ug/L	0.29429	0.1338 ug/L	0.29429	219.87%
QC value within limits for Cu 324.752		Recovery =	Not calculated			
Fe 238.204†	0.7	4.719 ug/L	5.4441	4.719 ug/L	5.4441	115.36%
QC value within limits for Fe 238.204		Recovery =	Not calculated			
K 766.490†	63.4	36.08 ug/L	12.580	36.08 ug/L	12.580	34.87%
QC value within limits for K 766.490		Recovery =	Not calculated			
Mg 285.213†	0.8	0.1089 ug/L	0.18918	0.1089 ug/L	0.18918	173.68%
QC value within limits for Mg 285.213		Recovery =	Not calculated			
Mn 257.610†	9.2	0.0628 ug/L	0.04181	0.0628 ug/L	0.04181	66.56%
QC value within limits for Mn 257.610		Recovery =	Not calculated			
Mo 202.031†	3.9	7.913 ug/L	0.5836	7.913 ug/L	0.5836	7.37%
QC value within limits for Mo 202.031		Recovery =	Not calculated			
Na 589.592†	21.4	3.746 ug/L	3.1989	3.746 ug/L	3.1989	85.39%
QC value within limits for Na 589.592		Recovery =	Not calculated			
Ni 231.604†	-1.5	-0.8094 ug/L	1.21970	-0.8094 ug/L	1.21970	150.69%
QC value within limits for Ni 231.604		Recovery =	Not calculated			
Pb 220.353†	0.3	0.7438 ug/L	11.45580	0.7438 ug/L	11.45580	>999.9%
QC value within limits for Pb 220.353		Recovery =	Not calculated			
Sb 206.836†	0.3	1.217 ug/L	5.9174	1.217 ug/L	5.9174	486.37%
QC value within limits for Sb 206.836		Recovery =	Not calculated			
Se 196.026†	2.4	28.17 ug/L	14.804	28.17 ug/L	14.804	52.55%
QC value within limits for Se 196.026		Recovery =	Not calculated			
SiO2 251.603†	1.0	0.2402 ug/L	0.40663	0.2402 ug/L	0.40663	169.28%
QC value within limits for SiO2 251.603		Recovery =	Not calculated			
Sr 421.552†	56.5	0.037 ug/L	0.0024	0.037 ug/L	0.0024	6.46%
QC value within limits for Sr 421.552		Recovery =	Not calculated			
Ti 334.940†	103.6	0.490 ug/L	0.0420	0.490 ug/L	0.0420	8.56%
QC value within limits for Ti 334.940		Recovery =	Not calculated			
Tl 190.801†	-0.2	-0.744 ug/L	10.3809	-0.744 ug/L	10.3809	>999.9%
QC value within limits for Tl 190.801		Recovery =	Not calculated			
V 290.880†	-43.6	-1.502 ug/L	0.3804	-1.502 ug/L	0.3804	25.32%
QC value within limits for V 290.880		Recovery =	Not calculated			

Zn 206.200†	1.7	1.545 ug/L	0.5561	1.545 ug/L	0.5561	35.99%
-------------	-----	------------	--------	------------	--------	--------

QC value within limits for Zn 206.200 Recovery = Not calculated

All analyte(s) passed QC.

Sequence No.: 5
 Sample ID: SEQ-CRL
 Analyst:
 Initial Sample Wt:
 Dilution:

Autosampler Location: 11
 Date Collected: 11/29/2010 7:50:20 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: SEQ-CRL

Analyte	Back Pressure	Flow
All	115.0 kPa	0.80 L/min

Mean Data: SEQ-CRL

Analyte	Mean Corrected Intensity	Calib Conc. Units	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2521949.6	102.9 %	0.80			0.78%
Sc Radial	313438.7	104.6 %	0.99			0.95%
Ag 328.068†	631.0	10.37 ug/L	0.229	10.37 ug/L	0.229	2.21%
QC value within limits for Ag 328.068			Recovery = 103.70%			
Al 396.153†	471.5	101.2 ug/L	1.57	101.2 ug/L	1.57	1.55%
QC value within limits for Al 396.153			Recovery = 101.21%			
As 193.696†	5.5	56.49 ug/L	14.366	56.49 ug/L	14.366	25.43%
QC value within limits for As 193.696			Recovery = 112.98%			
Ba 233.527†	125.5	10.43 ug/L	0.179	10.43 ug/L	0.179	1.71%
QC value within limits for Ba 233.527			Recovery = 104.31%			
Be 313.107†	1401.8	5.109 ug/L	0.0470	5.109 ug/L	0.0470	0.92%
QC value within limits for Be 313.107			Recovery = 102.18%			
B 249.677†	1744.4	267.9 ug/L	2.01	267.9 ug/L	2.01	0.75%
QC value within limits for B 249.677			Recovery = 107.15%			
Ca 317.933†	737.4	254.2 ug/L	4.23	254.2 ug/L	4.23	1.66%
QC value within limits for Ca 317.933			Recovery = 101.68%			
Cd 214.440†	36.5	11.15 ug/L	0.270	11.15 ug/L	0.270	2.42%
QC value within limits for Cd 214.440			Recovery = 111.51%			
Co 228.616†	45.6	10.76 ug/L	0.160	10.76 ug/L	0.160	1.49%
QC value within limits for Co 228.616			Recovery = 107.57%			
Cr 267.716†	60.6	9.655 ug/L	0.0966	9.655 ug/L	0.0966	1.00%
QC value within limits for Cr 267.716			Recovery = 96.55%			
Cu 324.752†	948.9	9.923 ug/L	0.3443	9.923 ug/L	0.3443	3.47%
QC value within limits for Cu 324.752			Recovery = 99.23%			
Fe 238.204†	14.1	99.55 ug/L	10.617	99.55 ug/L	10.617	10.67%
QC value within limits for Fe 238.204			Recovery = 99.55%			
K 766.490†	1849.5	1049 ug/L	19.7	1049 ug/L	19.7	1.88%
QC value within limits for K 766.490			Recovery = 104.86%			
Mg 285.213†	6684.7	1041 ug/L	10.3	1041 ug/L	10.3	0.99%
QC value within limits for Mg 285.213			Recovery = 104.11%			
Mn 257.610†	1420.2	10.33 ug/L	0.131	10.33 ug/L	0.131	1.27%
QC value within limits for Mn 257.610			Recovery = 103.25%			
Mo 202.031†	5.2	10.61 ug/L	1.032	10.61 ug/L	1.032	9.73%
QC value within limits for Mo 202.031			Recovery = 106.11%			
Na 589.592†	5679.0	1041 ug/L	15.9	1041 ug/L	15.9	1.53%
QC value within limits for Na 589.592			Recovery = 104.10%			
Ni 231.604†	19.5	10.84 ug/L	0.221	10.84 ug/L	0.221	2.04%
QC value within limits for Ni 231.604			Recovery = 108.44%			
Pb 220.353†	10.0	22.27 ug/L	5.098	22.27 ug/L	5.098	22.89%
QC value within limits for Pb 220.353			Recovery = 74.25%			
Sb 206.836†	10.0	38.44 ug/L	10.760	38.44 ug/L	10.760	27.99%
QC value within limits for Sb 206.836			Recovery = 76.88%			
Se 196.026†	10.7	126.3 ug/L	21.19	126.3 ug/L	21.19	16.77%
QC value within limits for Se 196.026			Recovery = 126.32%			
SiO2 251.603†	818.2	246.2 ug/L	1.90	246.2 ug/L	1.90	0.77%
QC value within limits for SiO2 251.603			Recovery = 98.47%			
Sr 421.552†	16151.2	10.58 ug/L	0.136	10.58 ug/L	0.136	1.28%
QC value within limits for Sr 421.552			Recovery = 105.80%			
Ti 334.940†	10816.7	51.17 ug/L	0.466	51.17 ug/L	0.466	0.91%
QC value within limits for Ti 334.940			Recovery = 102.33%			
Tl 190.801†	13.9	57.24 ug/L	10.961	57.24 ug/L	10.961	19.15%
QC value within limits for Tl 190.801			Recovery = 114.48%			
V 290.880†	1369.7	48.27 ug/L	0.398	48.27 ug/L	0.398	0.82%
QC value within limits for V 290.880			Recovery = 96.55%			

Zn 206.200†	56.6	53.21 ug/L	0.836	53.21 ug/L	0.836	1.57%
-------------	------	------------	-------	------------	-------	-------

QC value within limits for Zn 206.200 Recovery = 106.42%

All analyte(s) passed QC.

Sequence No.: 6
 Sample ID: SEQ-IFA
 Analyst:
 Initial Sample Wt:
 Dilution:

Autosampler Location: 12
 Date Collected: 11/29/2010 7:53:25 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: SEQ-IFA

Analyte	Back Pressure	Flow
All	116.0 kPa	0.80 L/min

Mean Data: SEQ-IFA

Analyte	Mean Corrected Intensity	Calib Conc. Units	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2347455.9	95.78 %	0.340			0.35%
Sc Radial	304131.6	101.5 %	0.46			0.45%
Ag 328.068†	-554.5	6.161 ug/L	0.3828	6.161 ug/L	0.3828	6.21%
QC value within limits for Ag 328.068		Recovery = Not calculated				
Al 396.153†	285649.3	61260 ug/L	387.6	61260 ug/L	387.6	0.63%
QC value within limits for Al 396.153		Recovery = 102.09%				
As 193.696†	-11.3	-5.916 ug/L	30.6957	-5.916 ug/L	30.6957	518.88%
QC value within limits for As 193.696		Recovery = Not calculated				
Ba 233.527†	29.2	-2.802 ug/L	0.5102	-2.802 ug/L	0.5102	18.21%
QC value within limits for Ba 233.527		Recovery = Not calculated				
Be 313.107†	-13.2	-0.7856 ug/L	0.05393	-0.7856 ug/L	0.05393	6.86%
QC value within limits for Be 313.107		Recovery = Not calculated				
B 249.677†	-832.6	-127.8 ug/L	1.07	-127.8 ug/L	1.07	0.84%
QC value within limits for B 249.677		Recovery = Not calculated				
Ca 317.933†	841846.7	292100 ug/L	1223.7	292100 ug/L	1223.7	0.42%
QC value within limits for Ca 317.933		Recovery = 97.36%				
Cd 214.440†	-12.6	-3.392 ug/L	0.4977	-3.392 ug/L	0.4977	14.68%
QC value within limits for Cd 214.440		Recovery = Not calculated				
Co 228.616†	31.2	4.065 ug/L	1.4661	4.065 ug/L	1.4661	36.07%
QC value within limits for Co 228.616		Recovery = Not calculated				
Cr 267.716†	-21.3	2.863 ug/L	0.4843	2.863 ug/L	0.4843	16.91%
QC value within limits for Cr 267.716		Recovery = Not calculated				
Cu 324.752†	-97.1	-0.8254 ug/L	0.09334	-0.8254 ug/L	0.09334	11.31%
QC value within limits for Cu 324.752		Recovery = Not calculated				
Fe 238.204†	32936.1	235100 ug/L	2766.2	235100 ug/L	2766.2	1.18%
QC value within limits for Fe 238.204		Recovery = 94.03%				
K 766.490†	89.2	-50.85 ug/L	20.768	-50.85 ug/L	20.768	40.84%
QC value within limits for K 766.490		Recovery = Not calculated				
Mg 285.213†	922531.8	143700 ug/L	767.3	143700 ug/L	767.3	0.53%
QC value within limits for Mg 285.213		Recovery = 95.83%				
Mn 257.610†	541.4	1.443 ug/L	0.2704	1.443 ug/L	0.2704	18.74%
QC value within limits for Mn 257.610		Recovery = Not calculated				
Mo 202.031†	2.3	5.341 ug/L	6.7121	5.341 ug/L	6.7121	125.68%
QC value within limits for Mo 202.031		Recovery = Not calculated				
Na 589.592†	278637.9	51070 ug/L	512.0	51070 ug/L	512.0	1.00%
QC value within limits for Na 589.592		Recovery = 102.14%				
Ni 231.604†	32.8	3.157 ug/L	2.2443	3.157 ug/L	2.2443	71.09%
QC value within limits for Ni 231.604		Recovery = Not calculated				
Pb 220.353†	-23.9	-18.91 ug/L	14.108	-18.91 ug/L	14.108	74.60%
QC value within limits for Pb 220.353		Recovery = Not calculated				
Sb 206.836†	16.4	55.46 ug/L	16.070	55.46 ug/L	16.070	28.98%
QC value within limits for Sb 206.836		Recovery = Not calculated				
Se 196.026†	2.7	68.66 ug/L	51.345	68.66 ug/L	51.345	74.78%
QC value within limits for Se 196.026		Recovery = Not calculated				
SiO2 251.603†	-58.4	62.95 ug/L	1.690	62.95 ug/L	1.690	2.68%
QC value within limits for SiO2 251.603		Recovery = Not calculated				
Sr 421.552†	4322.5	2.831 ug/L	0.0138	2.831 ug/L	0.0138	0.49%
QC value within limits for Sr 421.552		Recovery = Not calculated				
Ti 334.940†	42.9	0.203 ug/L	0.0673	0.203 ug/L	0.0673	33.14%
QC value within limits for Ti 334.940		Recovery = Not calculated				
Tl 190.801†	0.9	-7.823 ug/L	3.4949	-7.823 ug/L	3.4949	44.67%
QC value within limits for Tl 190.801		Recovery = Not calculated				
V 290.880†	704.3	-5.538 ug/L	1.6638	-5.538 ug/L	1.6638	30.05%
QC value within limits for V 290.880		Recovery = Not calculated				

Zn 206.200†	4.0	2.843 ug/L	2.3613	2.843 ug/L	2.3613	83.05%
-------------	-----	------------	--------	------------	--------	--------

QC value within limits for Zn 206.200 Recovery = Not calculated
All analyte(s) passed QC.

Sequence No.: 7
 Sample ID: SEQ-IFB
 Analyst:
 Initial Sample Wt:
 Dilution:

Autosampler Location: 13
 Date Collected: 11/29/2010 7:57:08 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: SEQ-IFB

Analyte	Back Pressure	Flow
All	115.0 kPa	0.80 L/min

Mean Data: SEQ-IFB

Analyte	Mean Corrected Intensity	Calib Conc. Units	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2332430.2	95.17 %	0.614			0.65%
Sc Radial	307310.4	102.5 %	1.08			1.05%
Ag 328.068†	18371.4	316.2 ug/L	1.89	316.2 ug/L	1.89	0.60%
QC value within limits for Ag 328.068 Recovery = 105.39%						
Al 396.153†	282069.6	60500 ug/L	1208.3	60500 ug/L	1208.3	2.00%
QC value within limits for Al 396.153 Recovery = 100.83%						
As 193.696†	85.1	969.9 ug/L	34.28	969.9 ug/L	34.28	3.53%
QC value within limits for As 193.696 Recovery = 96.99%						
Ba 233.527†	3541.2	289.4 ug/L	1.30	289.4 ug/L	1.30	0.45%
QC value within limits for Ba 233.527 Recovery = 96.46%						
Be 313.107†	26499.8	95.87 ug/L	0.684	95.87 ug/L	0.684	0.71%
QC value within limits for Be 313.107 Recovery = 95.87%						
B 249.677†	2396.4	368.0 ug/L	9.18	368.0 ug/L	9.18	2.50%
QC value less than the lower limit for B 249.677 Recovery = 73.60%						
Ca 317.933†	830194.1	288000 ug/L	6931.1	288000 ug/L	6931.1	2.41%
QC value within limits for Ca 317.933 Recovery = 96.00%						
Cd 214.440†	948.3	289.9 ug/L	0.72	289.9 ug/L	0.72	0.25%
QC value within limits for Cd 214.440 Recovery = 96.64%						
Co 228.616†	1261.1	293.6 ug/L	0.70	293.6 ug/L	0.70	0.24%
QC value within limits for Co 228.616 Recovery = 97.87%						
Cr 267.716†	1796.7	291.8 ug/L	1.85	291.8 ug/L	1.85	0.64%
QC value within limits for Cr 267.716 Recovery = 97.28%						
Cu 324.752†	29396.5	307.1 ug/L	3.91	307.1 ug/L	3.91	1.27%
QC value within limits for Cu 324.752 Recovery = 102.36%						
Fe 238.204†	32504.6	232000 ug/L	3386.6	232000 ug/L	3386.6	1.46%
QC value within limits for Fe 238.204 Recovery = 92.79%						
K 766.490†	35563.8	19980 ug/L	211.0	19980 ug/L	211.0	1.06%
QC value within limits for K 766.490 Recovery = 99.90%						
Mg 285.213†	913643.9	142400 ug/L	2913.5	142400 ug/L	2913.5	2.05%
QC value within limits for Mg 285.213 Recovery = 94.91%						
Mn 257.610†	27069.5	194.7 ug/L	1.70	194.7 ug/L	1.70	0.87%
QC value within limits for Mn 257.610 Recovery = 97.36%						
Mo 202.031†	144.1	293.9 ug/L	7.35	293.9 ug/L	7.35	2.50%
QC value within limits for Mo 202.031 Recovery = 97.98%						
Na 589.592†	275042.1	50360 ug/L	843.2	50360 ug/L	843.2	1.67%
QC value within limits for Na 589.592 Recovery = 100.72%						
Ni 231.604†	536.4	283.6 ug/L	2.41	283.6 ug/L	2.41	0.85%
QC value within limits for Ni 231.604 Recovery = 94.52%						
Pb 220.353†	411.9	959.4 ug/L	3.45	959.4 ug/L	3.45	0.36%
QC value within limits for Pb 220.353 Recovery = 95.94%						
Sb 206.836†	249.6	949.2 ug/L	23.63	949.2 ug/L	23.63	2.49%
QC value within limits for Sb 206.836 Recovery = 94.92%						
Se 196.026†	45.1	573.1 ug/L	80.83	573.1 ug/L	80.83	14.11%
QC value within limits for Se 196.026 Recovery = 114.61%						
SiO2 251.603†	1656.8	520.2 ug/L	4.37	520.2 ug/L	4.37	0.84%
QC value within limits for SiO2 251.603 Recovery = 104.04%						
Sr 421.552†	1552517.5	1017 ug/L	2.1	1017 ug/L	2.1	0.21%
QC value within limits for Sr 421.552 Recovery = 101.70%						
Ti 334.940†	216009.2	1022 ug/L	8.7	1022 ug/L	8.7	0.85%
QC value within limits for Ti 334.940 Recovery = 102.18%						
Tl 190.801†	243.7	982.5 ug/L	5.73	982.5 ug/L	5.73	0.58%
QC value within limits for Tl 190.801 Recovery = 98.25%						
V 290.880†	9257.8	295.3 ug/L	3.84	295.3 ug/L	3.84	1.30%
QC value within limits for V 290.880 Recovery = 98.43%						

Zn 206.200† 297.6 275.4 ug/L 0.84 275.4 ug/L 0.84 0.30%
QC value within limits for Zn 206.200 Recovery = 91.82%
QC Failed. Continue with analysis.

Sequence No.: 8
 Sample ID: SEQ-CCV
 Analyst:
 Initial Sample Wt:
 Dilution:

Autosampler Location: 3
 Date Collected: 11/29/2010 8:00:53 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: SEQ-CCV

Analyte	Back Pressure	Flow
All	116.0 kPa	0.80 L/min

Mean Data: SEQ-CCV

Analyte	Mean Corrected Intensity	Calib Conc. Units	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2498709.2	102.0 %	0.51			0.50%
Sc Radial	314554.4	104.9 %	0.64			0.61%
Ag 328.068†	15533.8	256.8 ug/L	2.25	256.8 ug/L	2.25	0.88%
QC value within limits for Ag 328.068		Recovery = 102.72%				
Al 396.153†	59311.2	12730 ug/L	40.3	12730 ug/L	40.3	0.32%
QC value within limits for Al 396.153		Recovery = 101.84%				
As 193.696†	250.5	2560 ug/L	36.8	2560 ug/L	36.8	1.44%
QC value within limits for As 193.696		Recovery = 102.42%				
Ba 233.527†	6033.8	501.4 ug/L	4.13	501.4 ug/L	4.13	0.82%
QC value within limits for Ba 233.527		Recovery = 100.27%				
Be 313.107†	141242.6	515.2 ug/L	0.84	515.2 ug/L	0.84	0.16%
QC value within limits for Be 313.107		Recovery = 103.03%				
B 249.677†	32778.8	5033 ug/L	38.5	5033 ug/L	38.5	0.77%
QC value within limits for B 249.677		Recovery = 100.67%				
Ca 317.933†	36761.4	12690 ug/L	46.8	12690 ug/L	46.8	0.37%
QC value within limits for Ca 317.933		Recovery = 101.55%				
Cd 214.440†	1692.8	516.8 ug/L	5.91	516.8 ug/L	5.91	1.14%
QC value within limits for Cd 214.440		Recovery = 103.35%				
Co 228.616†	2193.5	516.2 ug/L	3.53	516.2 ug/L	3.53	0.68%
QC value within limits for Co 228.616		Recovery = 103.25%				
Cr 267.716†	15635.2	2479 ug/L	35.0	2479 ug/L	35.0	1.41%
QC value within limits for Cr 267.716		Recovery = 99.16%				
Cu 324.752†	98136.5	1024 ug/L	2.7	1024 ug/L	2.7	0.26%
QC value within limits for Cu 324.752		Recovery = 102.37%				
Fe 238.204†	1773.5	12630 ug/L	46.4	12630 ug/L	46.4	0.37%
QC value within limits for Fe 238.204		Recovery = 101.05%				
K 766.490†	44718.6	25270 ug/L	99.3	25270 ug/L	99.3	0.39%
QC value within limits for K 766.490		Recovery = 101.07%				
Mg 285.213†	82321.9	12810 ug/L	48.6	12810 ug/L	48.6	0.38%
QC value within limits for Mg 285.213		Recovery = 102.51%				
Mn 257.610†	139203.4	1015 ug/L	8.3	1015 ug/L	8.3	0.82%
QC value within limits for Mn 257.610		Recovery = 101.51%				
Mo 202.031†	247.5	501.7 ug/L	0.92	501.7 ug/L	0.92	0.18%
QC value within limits for Mo 202.031		Recovery = 100.35%				
Na 589.592†	69751.6	12700 ug/L	44.9	12700 ug/L	44.9	0.35%
QC value within limits for Na 589.592		Recovery = 101.61%				
Ni 231.604†	4602.0	2551 ug/L	35.2	2551 ug/L	35.2	1.38%
QC value within limits for Ni 231.604		Recovery = 102.03%				
Pb 220.353†	1148.3	2586 ug/L	27.9	2586 ug/L	27.9	1.08%
QC value within limits for Pb 220.353		Recovery = 103.44%				
Sb 206.836†	655.8	2494 ug/L	12.5	2494 ug/L	12.5	0.50%
QC value within limits for Sb 206.836		Recovery = 99.76%				
Se 196.026†	225.8	2665 ug/L	46.4	2665 ug/L	46.4	1.74%
QC value within limits for Se 196.026		Recovery = 106.59%				
SiO2 251.603†	33482.7	10090 ug/L	99.5	10090 ug/L	99.5	0.99%
QC value within limits for SiO2 251.603		Recovery = 100.87%				
Sr 421.552†	794593.6	520.5 ug/L	0.43	520.5 ug/L	0.43	0.08%
QC value within limits for Sr 421.552		Recovery = 104.10%				
Ti 334.940†	107257.5	507.4 ug/L	0.36	507.4 ug/L	0.36	0.07%
QC value within limits for Ti 334.940		Recovery = 101.47%				
Tl 190.801†	634.3	2594 ug/L	1.8	2594 ug/L	1.8	0.07%
QC value within limits for Tl 190.801		Recovery = 103.75%				
V 290.880†	28318.3	998.4 ug/L	11.63	998.4 ug/L	11.63	1.16%
QC value within limits for V 290.880		Recovery = 99.84%				

Zn 206.200†	2716.3	2561 ug/L	26.9	2561 ug/L	26.9	1.05%
-------------	--------	-----------	------	-----------	------	-------

QC value within limits for Zn 206.200 Recovery = 102.42%

All analyte(s) passed QC.

Sequence No.: 9
 Sample ID: SEQ-CCB
 Analyst:
 Initial Sample Wt:
 Dilution:

Autosampler Location: 1
 Date Collected: 11/29/2010 8:04:00 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: SEQ-CCB

Analyte	Back Pressure	Flow
All	116.0 kPa	0.80 L/min

Mean Data: SEQ-CCB

Analyte	Mean Corrected Intensity	Conc. Units	Calib	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2520617.8	102.8 %		0.61			0.59%
Sc Radial	311314.5	103.9 %		0.87			0.84%
Ag 328.068†	10.4	0.1771 ug/L		0.04272	0.1771 ug/L	0.04272	24.13%
QC value within limits for Ag 328.068			Recovery =	Not calculated			
Al 396.153†	5.2	0.8182 ug/L		0.67406	0.8182 ug/L	0.67406	82.38%
QC value within limits for Al 396.153			Recovery =	Not calculated			
As 193.696†	0.1	1.370 ug/L		9.0247	1.370 ug/L	9.0247	658.61%
QC value within limits for As 193.696			Recovery =	Not calculated			
Ba 233.527†	1.6	0.1150 ug/L		0.06847	0.1150 ug/L	0.06847	59.55%
QC value within limits for Ba 233.527			Recovery =	Not calculated			
Be 313.107†	6.8	0.0238 ug/L		0.05720	0.0238 ug/L	0.05720	240.21%
QC value within limits for Be 313.107			Recovery =	Not calculated			
B 249.677†	317.8	48.79 ug/L		4.184	48.79 ug/L	4.184	8.57%
QC value within limits for B 249.677			Recovery =	Not calculated			
Ca 317.933†	3.1	0.9202 ug/L		0.43621	0.9202 ug/L	0.43621	47.40%
QC value within limits for Ca 317.933			Recovery =	Not calculated			
Cd 214.440†	0.1	0.0193 ug/L		0.46680	0.0193 ug/L	0.46680	>999.9%
QC value within limits for Cd 214.440			Recovery =	Not calculated			
Co 228.616†	0.8	0.1820 ug/L		0.06618	0.1820 ug/L	0.06618	36.36%
QC value within limits for Co 228.616			Recovery =	Not calculated			
Cr 267.716†	-0.2	-0.0259 ug/L		0.30632	-0.0259 ug/L	0.30632	>999.9%
QC value within limits for Cr 267.716			Recovery =	Not calculated			
Cu 324.752†	57.9	0.6079 ug/L		0.04883	0.6079 ug/L	0.04883	8.03%
QC value within limits for Cu 324.752			Recovery =	Not calculated			
Fe 238.204†	5.4	38.73 ug/L		17.301	38.73 ug/L	17.301	44.67%
QC value within limits for Fe 238.204			Recovery =	Not calculated			
K 766.490†	59.3	33.35 ug/L		20.044	33.35 ug/L	20.044	60.10%
QC value within limits for K 766.490			Recovery =	Not calculated			
Mg 285.213†	9.9	1.510 ug/L		0.6200	1.510 ug/L	0.6200	41.07%
QC value within limits for Mg 285.213			Recovery =	Not calculated			
Mn 257.610†	9.9	0.0670 ug/L		0.01514	0.0670 ug/L	0.01514	22.62%
QC value within limits for Mn 257.610			Recovery =	Not calculated			
Mo 202.031†	3.4	6.930 ug/L		1.0673	6.930 ug/L	1.0673	15.40%
QC value within limits for Mo 202.031			Recovery =	Not calculated			
Na 589.592†	-14.8	-2.976 ug/L		4.7126	-2.976 ug/L	4.7126	158.38%
QC value within limits for Na 589.592			Recovery =	Not calculated			
Ni 231.604†	-0.0	-0.0097 ug/L		0.67369	-0.0097 ug/L	0.67369	>999.9%
QC value within limits for Ni 231.604			Recovery =	Not calculated			
Pb 220.353†	-2.3	-5.222 ug/L		11.8860	-5.222 ug/L	11.8860	227.59%
QC value within limits for Pb 220.353			Recovery =	Not calculated			
Sb 206.836†	6.5	25.14 ug/L		10.640	25.14 ug/L	10.640	42.32%
QC value within limits for Sb 206.836			Recovery =	Not calculated			
Se 196.026†	2.3	26.89 ug/L		46.334	26.89 ug/L	46.334	172.31%
QC value within limits for Se 196.026			Recovery =	Not calculated			
SiO2 251.603†	9.1	2.771 ug/L		0.1315	2.771 ug/L	0.1315	4.75%
QC value within limits for SiO2 251.603			Recovery =	Not calculated			
Sr 421.552†	56.6	0.037 ug/L		0.0064	0.037 ug/L	0.0064	17.32%
QC value within limits for Sr 421.552			Recovery =	Not calculated			
Ti 334.940†	60.8	0.288 ug/L		0.0346	0.288 ug/L	0.0346	12.03%
QC value within limits for Ti 334.940			Recovery =	Not calculated			
Tl 190.801†	-0.3	-1.506 ug/L		13.6811	-1.506 ug/L	13.6811	908.51%
QC value within limits for Tl 190.801			Recovery =	Not calculated			
V 290.880†	-50.4	-1.757 ug/L		0.4178	-1.757 ug/L	0.4178	23.78%
QC value within limits for V 290.880			Recovery =	Not calculated			

Zn 206.200†	1.9	1.741 ug/L	1.1960	1.741 ug/L	1.1960	68.70%
-------------	-----	------------	--------	------------	--------	--------

QC value within limits for Zn 206.200 Recovery = Not calculated

All analyte(s) passed QC.

Sequence No.: 10
 Sample ID: 1011109-BLK1
 Analyst: Walker
 Initial Sample Wt:
 Dilution:

Autosampler Location: 26
 Date Collected: 11/29/2010 8:07:03 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: 1011109-BLK1

Analyte Back Pressure Flow
 All 115.0 kPa 0.80 L/min

Mean Data: 1011109-BLK1

Analyte	Mean Corrected Intensity	Calib Conc. Units	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2514925.7	102.6 %	0.29			0.28%
Sc Radial	310742.0	103.7 %	1.26			1.22%
Ag 328.068†	27.7	0.4539 ug/L	0.24330	0.4539 ug/L	0.24330	53.60%
Al 396.153†	-1.9	-0.5729 ug/L	2.94080	-0.5729 ug/L	2.94080	513.29%
As 193.696†	0.4	4.055 ug/L	20.1992	4.055 ug/L	20.1992	498.11%
Ba 233.527†	1.3	0.1046 ug/L	0.05705	0.1046 ug/L	0.05705	54.55%
Be 313.107†	-10.4	-0.0384 ug/L	0.02084	-0.0384 ug/L	0.02084	54.24%
B 249.677†	150.7	23.15 ug/L	1.162	23.15 ug/L	1.162	5.02%
Ca 317.933†	-1.7	-0.6523 ug/L	3.98143	-0.6523 ug/L	3.98143	610.35%
Cd 214.440†	1.9	0.5848 ug/L	0.28358	0.5848 ug/L	0.28358	48.49%
Co 228.616†	2.6	0.6165 ug/L	0.93344	0.6165 ug/L	0.93344	151.41%
Cr 267.716†	-1.7	-0.2707 ug/L	0.24419	-0.2707 ug/L	0.24419	90.22%
Cu 324.752†	51.8	0.5410 ug/L	0.01414	0.5410 ug/L	0.01414	2.61%
Fe 238.204†	3.0	21.14 ug/L	8.557	21.14 ug/L	8.557	40.48%
K 766.490†	49.4	28.11 ug/L	7.328	28.11 ug/L	7.328	26.07%
Mg 285.213†	-4.5	-0.7155 ug/L	0.13602	-0.7155 ug/L	0.13602	19.01%
Mn 257.610†	-4.8	-0.0366 ug/L	0.01976	-0.0366 ug/L	0.01976	53.97%
Mo 202.031†	1.5	3.019 ug/L	2.1263	3.019 ug/L	2.1263	70.42%
Na 589.592†	-48.7	-9.041 ug/L	3.3315	-9.041 ug/L	3.3315	36.85%
Ni 231.604†	0.9	0.4822 ug/L	1.62380	0.4822 ug/L	1.62380	336.75%
Pb 220.353†	-2.3	-5.252 ug/L	1.9795	-5.252 ug/L	1.9795	37.69%
Sb 206.836†	2.1	8.160 ug/L	2.8256	8.160 ug/L	2.8256	34.63%
Se 196.026†	0.7	8.503 ug/L	12.9807	8.503 ug/L	12.9807	152.66%
SiO2 251.603†	128.3	38.65 ug/L	0.447	38.65 ug/L	0.447	1.16%
Sr 421.552†	-7.5	-0.005 ug/L	0.0063	-0.005 ug/L	0.0063	128.69%
Ti 334.940†	30.0	0.142 ug/L	0.0074	0.142 ug/L	0.0074	5.19%
Tl 190.801†	-0.0	-0.105 ug/L	8.2143	-0.105 ug/L	8.2143	>999.9%
V 290.880†	-55.7	-1.957 ug/L	0.1429	-1.957 ug/L	0.1429	7.30%
Zn 206.200†	1.0	0.909 ug/L	0.3356	0.909 ug/L	0.3356	36.92%

Sequence No.: 11
 Sample ID: 1011109-BS1
 Analyst: Walker
 Initial Sample Wt:
 Dilution:

Autosampler Location: 27
 Date Collected: 11/29/2010 8:10:06 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: 1011109-BS1

Analyte Back Pressure Flow
 All 117.0 kPa 0.80 L/min

Mean Data: 1011109-BS1

Analyte	Mean Corrected Intensity	Conc. Units	Calib	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2514144.0	102.6 %		0.56			0.54%
Sc Radial	314996.6	105.1 %		0.24			0.23%
Ag 328.068†	6120.0	101.3 ug/L		0.64	101.3 ug/L	0.64	0.63%
Al 396.153†	48787.9	10480 ug/L		85.9	10480 ug/L	85.9	0.82%
As 193.696†	7.5	87.86 ug/L		27.124	87.86 ug/L	27.124	30.87%
Ba 233.527†	1195.3	98.91 ug/L		0.589	98.91 ug/L	0.589	0.60%
Be 313.107†	27795.5	101.4 ug/L		0.60	101.4 ug/L	0.60	0.59%
B 249.677†	0.8	0.1187 ug/L		0.57101	0.1187 ug/L	0.57101	481.13%
Ca 317.933†	29453.3	10200 ug/L		25.0	10200 ug/L	25.0	0.25%
Cd 214.440†	336.6	102.8 ug/L		0.79	102.8 ug/L	0.79	0.77%
Co 228.616†	422.6	99.43 ug/L		1.009	99.43 ug/L	1.009	1.01%
Cr 267.716†	614.3	98.04 ug/L		1.234	98.04 ug/L	1.234	1.26%
Cu 324.752†	9277.6	96.99 ug/L		0.513	96.99 ug/L	0.513	0.53%
Fe 238.204†	1450.7	10340 ug/L		21.1	10340 ug/L	21.1	0.20%
K 766.490†	18360.9	10380 ug/L		50.6	10380 ug/L	50.6	0.49%
Mg 285.213†	66941.9	10430 ug/L		82.0	10430 ug/L	82.0	0.79%
Mn 257.610†	13720.8	99.77 ug/L		0.243	99.77 ug/L	0.243	0.24%
Mo 202.031†	46.2	94.23 ug/L		6.446	94.23 ug/L	6.446	6.84%
Na 589.592†	57310.4	10500 ug/L		55.1	10500 ug/L	55.1	0.52%
Ni 231.604†	176.7	97.75 ug/L		0.302	97.75 ug/L	0.302	0.31%
Pb 220.353†	40.5	90.57 ug/L		1.963	90.57 ug/L	1.963	2.17%
Sb 206.836†	19.4	71.78 ug/L		3.558	71.78 ug/L	3.558	4.96%
Se 196.026†	47.1	557.2 ug/L		39.76	557.2 ug/L	39.76	7.14%
SiO2 251.603†	49.5	-3.185 ug/L		0.7157	-3.185 ug/L	0.7157	22.48%
Sr 421.552†	819879.6	537.1 ug/L		0.78	537.1 ug/L	0.78	0.14%
Ti 334.940†	-43.0	-0.203 ug/L		0.0096	-0.203 ug/L	0.0096	4.70%
Tl 190.801†	25.9	103.7 ug/L		2.14	103.7 ug/L	2.14	2.06%
V 290.880†	2704.1	92.95 ug/L		0.288	92.95 ug/L	0.288	0.31%
Zn 206.200†	108.4	100.4 ug/L		1.53	100.4 ug/L	1.53	1.53%

Matrix Recovery Check: 1011109-BS1

Analyte	Expected Conc.	Measured Conc.	Std. Dev.	Units	Recovery (%)
Al 396.153	10100	10480	85.917	ug/L	103.8
Ca 317.933	10100	10200	24.999	ug/L	101.0
Fe 238.204	10120	10340	21.085	ug/L	102.2
K 766.490	10130	10380	50.559	ug/L	102.4
Mg 285.213	10100	10430	82.026	ug/L	103.3
Na 589.592	10090	10500	55.061	ug/L	104.1
Ag 328.068	100.5	101.3	0.642	ug/L	100.9
As 193.696	104.1	87.86	27.124	ug/L	83.8
Ba 233.527	100.1	98.91	0.589	ug/L	98.8
Be 313.107	99.96	101.4	0.596	ug/L	101.4
Cd 214.440	100.6	102.8	0.788	ug/L	102.2
Co 228.616	100.6	99.43	1.009	ug/L	98.8
Cr 267.716	99.73	98.04	1.234	ug/L	98.3
Cu 324.752	100.5	96.99	0.513	ug/L	96.4
Mn 257.610	99.96	99.77	0.243	ug/L	99.8
Mo 202.031	103.0	94.23	6.446	ug/L	91.2
Ni 231.604	100.5	97.75	0.302	ug/L	97.3
Pb 220.353	94.75	90.57	1.963	ug/L	95.8
Sb 206.836	108.2	71.78	3.558	ug/L	63.6
Se 196.026	508.5	557.2	39.761	ug/L	109.7

SiO2 251.603	5039	-3.185	0.716	ug/L	-0.8
Sr 421.552	500.0	537.1	0.778	ug/L	107.4
Tl 190.801	99.90	103.7	2.141	ug/L	103.8
V 290.880	98.04	92.95	0.288	ug/L	94.9
Zn 206.200	100.9	100.4	1.534	ug/L	99.5

Sequence No.: 12
 Sample ID: C101104-03
 Analyst: Walker
 Initial Sample Wt:
 Dilution:

Autosampler Location: 28
 Date Collected: 11/29/2010 8:13:10 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: C101104-03

Analyte Back Pressure Flow
 All 115.0 kPa 0.80 L/min

Mean Data: C101104-03

Analyte	Mean Corrected Intensity	Conc. Units	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2488888.7	101.6 %	0.34			0.34%
Sc Radial	310758.2	103.7 %	1.00			0.96%
Ag 328.068†	3.8	0.7664 ug/L	0.14018	0.7664 ug/L	0.14018	18.29%
Al 396.153†	108.8	15.53 ug/L	0.909	15.53 ug/L	0.909	5.85%
As 193.696†	-7.3	-74.41 ug/L	7.856	-74.41 ug/L	7.856	10.56%
Ba 233.527†	812.4	67.18 ug/L	0.378	67.18 ug/L	0.378	0.56%
Be 313.107†	-5.7	-0.1036 ug/L	0.00919	-0.1036 ug/L	0.00919	8.88%
B 249.677†	-27.9	-4.283 ug/L	0.0612	-4.283 ug/L	0.0612	1.43%
Ca 317.933†	122421.0	42470 ug/L	255.6	42470 ug/L	255.6	0.60%
Cd 214.440†	0.6	0.2142 ug/L	0.32497	0.2142 ug/L	0.32497	151.72%
Co 228.616†	1.6	0.5661 ug/L	0.25157	0.5661 ug/L	0.25157	44.44%
Cr 267.716†	-2.9	0.1090 ug/L	0.11235	0.1090 ug/L	0.11235	103.10%
Cu 324.752†	89.6	1.109 ug/L	0.1740	1.109 ug/L	0.1740	15.68%
Fe 238.204†	5.4	32.72 ug/L	23.982	32.72 ug/L	23.982	73.29%
K 766.490†	1163.9	650.7 ug/L	8.15	650.7 ug/L	8.15	1.25%
Mg 285.213†	41318.1	6435 ug/L	30.5	6435 ug/L	30.5	0.47%
Mn 257.610†	1998.5	14.32 ug/L	0.079	14.32 ug/L	0.079	0.55%
Mo 202.031†	4.1	7.834 ug/L	3.5939	7.834 ug/L	3.5939	45.88%
Na 589.592†	14923.6	2733 ug/L	25.8	2733 ug/L	25.8	0.94%
Ni 231.604†	-0.3	-0.2071 ug/L	1.27420	-0.2071 ug/L	1.27420	615.18%
Pb 220.353†	-5.3	-12.01 ug/L	11.552	-12.01 ug/L	11.552	96.22%
Sb 206.836†	-0.2	-3.298 ug/L	10.5790	-3.298 ug/L	10.5790	320.78%
Se 196.026†	10.1	116.9 ug/L	13.28	116.9 ug/L	13.28	11.36%
SiO2 251.603†	19350.3	5836 ug/L	56.1	5836 ug/L	56.1	0.96%
Sr 421.552†	563087.0	368.8 ug/L	0.38	368.8 ug/L	0.38	0.10%
Ti 334.940†	-26.9	-0.127 ug/L	0.0360	-0.127 ug/L	0.0360	28.24%
Tl 190.801†	-0.6	-5.384 ug/L	7.1372	-5.384 ug/L	7.1372	132.55%
V 290.880†	-40.4	-2.801 ug/L	0.4172	-2.801 ug/L	0.4172	14.90%
Zn 206.200†	3.4	2.661 ug/L	1.4713	2.661 ug/L	1.4713	55.29%

Sequence No.: 13
 Sample ID: 1011109-DUP1
 Analyst: Walker
 Initial Sample Wt:
 Dilution:

Autosampler Location: 29
 Date Collected: 11/29/2010 8:16:14 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: 1011109-DUP1

Analyte Back Pressure Flow
 All 115.0 kPa 0.80 L/min

Mean Data: 1011109-DUP1

Analyte	Mean Corrected Intensity	Conc. Units	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2509161.0	102.4 %	0.32			0.32%
Sc Radial	310833.0	103.7 %	0.96			0.93%
Ag 328.068†	14.4	0.9276 ug/L	0.31755	0.9276 ug/L	0.31755	34.23%
Al 396.153†	103.2	14.41 ug/L	1.689	14.41 ug/L	1.689	11.72%
As 193.696†	-4.6	-47.27 ug/L	30.165	-47.27 ug/L	30.165	63.81%
Ba 233.527†	810.1	67.00 ug/L	0.496	67.00 ug/L	0.496	0.74%
Be 313.107†	-15.4	-0.1385 ug/L	0.06853	-0.1385 ug/L	0.06853	49.46%
B 249.677†	-40.3	-6.196 ug/L	0.4245	-6.196 ug/L	0.4245	6.85%
Ca 317.933†	121868.4	42280 ug/L	85.1	42280 ug/L	85.1	0.20%
Cd 214.440†	-0.0	0.0306 ug/L	0.77997	0.0306 ug/L	0.77997	>999.9%
Co 228.616†	3.3	0.9617 ug/L	0.61341	0.9617 ug/L	0.61341	63.79%
Cr 267.716†	-3.0	0.0994 ug/L	0.18068	0.0994 ug/L	0.18068	181.84%
Cu 324.752†	87.2	1.075 ug/L	0.1137	1.075 ug/L	0.1137	10.57%
Fe 238.204†	3.5	19.19 ug/L	22.975	19.19 ug/L	22.975	119.72%
K 766.490†	1146.0	640.5 ug/L	10.13	640.5 ug/L	10.13	1.58%
Mg 285.213†	41160.2	6410 ug/L	19.9	6410 ug/L	19.9	0.31%
Mn 257.610†	1978.0	14.18 ug/L	0.061	14.18 ug/L	0.061	0.43%
Mo 202.031†	3.4	6.521 ug/L	0.8372	6.521 ug/L	0.8372	12.84%
Na 589.592†	14909.2	2731 ug/L	18.3	2731 ug/L	18.3	0.67%
Ni 231.604†	-1.4	-0.8160 ug/L	1.21712	-0.8160 ug/L	1.21712	149.16%
Pb 220.353†	-6.0	-13.47 ug/L	11.002	-13.47 ug/L	11.002	81.70%
Sb 206.836†	-1.9	-9.656 ug/L	2.3892	-9.656 ug/L	2.3892	24.74%
Se 196.026†	6.7	76.03 ug/L	40.232	76.03 ug/L	40.232	52.92%
SiO2 251.603†	19219.2	5796 ug/L	16.9	5796 ug/L	16.9	0.29%
Sr 421.552†	563971.6	369.4 ug/L	0.52	369.4 ug/L	0.52	0.14%
Ti 334.940†	-24.9	-0.118 ug/L	0.0225	-0.118 ug/L	0.0225	19.11%
Tl 190.801†	-2.0	-10.87 ug/L	7.231	-10.87 ug/L	7.231	66.56%
V 290.880†	-66.7	-3.716 ug/L	0.0705	-3.716 ug/L	0.0705	1.90%
Zn 206.200†	3.6	2.839 ug/L	0.8608	2.839 ug/L	0.8608	30.32%

Duplicate Check: 1011109-DUP1

Analyte	Expected Conc.	Measured Conc.	Std. Dev.	Units	Difference (%)
Sc Radial	103.7	103.7	0.964	%	0.0
Al 396.153	15.53	14.41	1.689	ug/L	7.5
Ca 317.933	42470	42280	85.144	ug/L	0.5
Fe 238.204	32.72	19.19	22.975	ug/L	52.1
K 766.490	650.7	640.5	10.129	ug/L	1.6
Mg 285.213	6435	6410	19.896	ug/L	0.4
Na 589.592	2733	2731	18.274	ug/L	0.1
Sc Axial	101.6	102.4	0.325	%	0.8
Ag 328.068	0.7664	0.9276	0.318	ug/L	19.0
As 193.696	-74.41	-47.27	30.165	ug/L	-44.6
Ba 233.527	67.18	67.00	0.496	ug/L	0.3
Be 313.107	-0.1036	-0.1385	0.069	ug/L	-28.9
B 249.677	-4.283	-6.196	0.425	ug/L	-36.5
Cd 214.440	0.2142	0.0306	0.780	ug/L	150.0
Co 228.616	0.5661	0.9617	0.613	ug/L	51.8
Cr 267.716	0.1090	0.0994	0.181	ug/L	9.2
Cu 324.752	1.109	1.075	0.114	ug/L	3.1
Mn 257.610	14.32	14.18	0.061	ug/L	1.0
Mo 202.031	7.834	6.521	0.837	ug/L	18.3
Ni 231.604	-0.2071	-0.8160	1.217	ug/L	-119.0

Pb 220.353	-12.01	-13.47	11.002	ug/L	-11.5
Sb 206.836	-3.298	-9.656	2.389	ug/L	-98.2
Se 196.026	116.9	76.03	40.232	ug/L	42.4
SiO2 251.603	5836	5796	16.930	ug/L	0.7
Sr 421.552	368.8	369.4	0.520	ug/L	0.2
Ti 334.940	-0.127	-0.118	0.023	ug/L	-7.8
Tl 190.801	-5.384	-10.87	7.231	ug/L	-67.5
V 290.880	-2.801	-3.716	0.071	ug/L	-28.1
Zn 206.200	2.661	2.839	0.861	ug/L	6.5

Sequence No.: 14
 Sample ID: SEQ-SRD1 @5X
 Analyst: Walker
 Initial Sample Wt:
 Dilution: 5X

Autosampler Location: 30
 Date Collected: 11/29/2010 8:19:18 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: SEQ-SRD1 @5X

Analyte Back Pressure Flow
 All 115.0 kPa 0.80 L/min

Mean Data: SEQ-SRD1 @5X

Analyte	Mean Corrected Intensity	Conc. Units	Calib Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2514289.1	102.6 %	0.61			0.59%
Sc Radial	311578.4	104.0 %	0.67			0.65%
Ag 328.068†	24.2	0.5377 ug/L	0.09838	2.689 ug/L	0.4919	18.29%
Al 396.153†	11.6	0.8890 ug/L	1.78505	4.445 ug/L	8.9253	200.78%
As 193.696†	-1.1	-10.98 ug/L	3.451	-54.92 ug/L	17.255	31.42%
Ba 233.527†	163.7	13.53 ug/L	0.116	67.66 ug/L	0.582	0.86%
Be 313.107†	-9.2	-0.0503 ug/L	0.00779	-0.2513 ug/L	0.03896	15.50%
B 249.677†	1.8	0.2803 ug/L	0.30784	1.402 ug/L	1.5392	109.81%
Ca 317.933†	24288.2	8427 ug/L	78.8	42130 ug/L	393.8	0.93%
Cd 214.440†	0.1	0.0411 ug/L	0.18686	0.2056 ug/L	0.93432	454.38%
Co 228.616†	-1.0	-0.2080 ug/L	0.72339	-1.040 ug/L	3.6169	347.74%
Cr 267.716†	-1.8	-0.1774 ug/L	0.34162	-0.8871 ug/L	1.70811	192.56%
Cu 324.752†	40.5	0.4601 ug/L	0.18710	2.301 ug/L	0.9355	40.66%
Fe 238.204†	1.4	8.768 ug/L	4.4495	43.84 ug/L	22.248	50.75%
K 766.490†	244.2	136.2 ug/L	16.61	681.2 ug/L	83.04	12.19%
Mg 285.213†	8198.1	1277 ug/L	15.4	6384 ug/L	76.9	1.21%
Mn 257.610†	391.9	2.806 ug/L	0.0985	14.03 ug/L	0.492	3.51%
Mo 202.031†	1.7	3.446 ug/L	1.7987	17.23 ug/L	8.994	52.20%
Na 589.592†	2937.5	537.8 ug/L	6.17	2689 ug/L	30.9	1.15%
Ni 231.604†	0.5	0.2543 ug/L	0.91600	1.271 ug/L	4.5800	360.22%
Pb 220.353†	-0.2	-0.5081 ug/L	1.87089	-2.540 ug/L	9.3545	368.24%
Sb 206.836†	1.9	6.874 ug/L	5.4657	34.37 ug/L	27.328	79.51%
Se 196.026†	3.3	37.97 ug/L	44.671	189.9 ug/L	223.35	117.64%
SiO2 251.603†	3751.8	1131 ug/L	7.6	5657 ug/L	37.8	0.67%
Sr 421.552†	113553.5	74.38 ug/L	0.044	371.9 ug/L	0.22	0.06%
Ti 334.940†	2.4	0.011 ug/L	0.0106	0.056 ug/L	0.0530	94.24%
Tl 190.801†	2.0	7.599 ug/L	4.4596	37.99 ug/L	22.298	58.69%
V 290.880†	-85.1	-3.278 ug/L	0.3821	-16.39 ug/L	1.910	11.66%
Zn 206.200†	2.4	2.095 ug/L	0.6203	10.47 ug/L	3.102	29.62%

Dilution Check: SEQ-SRD1 @5X

Analyte	Expected Conc.	Measured Conc.	Std. Dev.	Units	Difference (%)
Sc Radial	20.74	104.0	0.672	%	401.3
Al 396.153	3.106	0.8890	1.785	ug/L	71.4
Ca 317.933	8495	8427	78.759	ug/L	0.8
Fe 238.204	6.545	8.768	4.450	ug/L	34.0
K 766.490	130.1	136.2	16.607	ug/L	4.7
Mg 285.213	1287	1277	15.387	ug/L	0.8
Na 589.592	546.7	537.8	6.170	ug/L	1.6
Sc Axial	20.31	102.6	0.607	%	405.1
Ag 328.068	0.1533	0.5377	0.098	ug/L	250.8
As 193.696	-14.88	-10.98	3.451	ug/L	-26.2
Ba 233.527	13.44	13.53	0.116	ug/L	0.7
Be 313.107	-0.0207	-0.0503	0.008	ug/L	-142.7
B 249.677	-0.8566	0.2803	0.308	ug/L	-132.7
Cd 214.440	0.0428	0.0411	0.187	ug/L	4.0
Co 228.616	0.1132	-0.2080	0.723	ug/L	283.7
Cr 267.716	0.0218	-0.1774	0.342	ug/L	914.0
Cu 324.752	0.2219	0.4601	0.187	ug/L	107.4
Mn 257.610	2.865	2.806	0.098	ug/L	2.0
Mo 202.031	1.567	3.446	1.799	ug/L	119.9
Ni 231.604	-0.0414	0.2543	0.916	ug/L	-713.8

Pb 220.353	-2.401	-0.5081	1.871	ug/L	-78.8
Sb 206.836	-0.6596	6.874	5.466	ug/L	-1142.2
Se 196.026	23.38	37.97	44.671	ug/L	62.4
SiO2 251.603	1167	1131	7.560	ug/L	3.1
Sr 421.552	73.77	74.38	0.044	ug/L	0.8
Ti 334.940	-0.025	0.011	0.011	ug/L	-144.2
Tl 190.801	-1.077	7.599	4.460	ug/L	-805.6
V 290.880	-0.560	-3.278	0.382	ug/L	-485.3
Zn 206.200	0.532	2.095	0.620	ug/L	293.5

Sequence No.: 15
 Sample ID: 1011109-MS1
 Analyst: Walker
 Initial Sample Wt:
 Dilution:

Autosampler Location: 31
 Date Collected: 11/29/2010 8:22:22 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: 1011109-MS1

Analyte	Back Pressure	Flow
All	115.0 kPa	0.80 L/min

Mean Data: 1011109-MS1

Analyte	Mean Corrected Intensity	Conc. Units	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2481674.5	101.3 %	0.59			0.59%
Sc Radial	309411.5	103.2 %	0.58			0.56%
Ag 328.068†	6288.2	104.8 ug/L	0.64	104.8 ug/L	0.64	0.61%
Al 396.153†	50346.0	10810 ug/L	98.9	10810 ug/L	98.9	0.92%
As 193.696†	2.9	40.03 ug/L	52.990	40.03 ug/L	52.990	132.39%
Ba 233.527†	2007.0	166.0 ug/L	1.18	166.0 ug/L	1.18	0.71%
Be 313.107†	29185.8	106.3 ug/L	0.82	106.3 ug/L	0.82	0.77%
B 249.677†	-114.1	-17.52 ug/L	0.309	-17.52 ug/L	0.309	1.77%
Ca 317.933†	149813.4	51960 ug/L	418.5	51960 ug/L	418.5	0.81%
Cd 214.440†	344.6	105.3 ug/L	0.45	105.3 ug/L	0.45	0.43%
Co 228.616†	430.4	101.4 ug/L	1.28	101.4 ug/L	1.28	1.26%
Cr 267.716†	622.5	99.91 ug/L	1.066	99.91 ug/L	1.066	1.07%
Cu 324.752†	9651.6	101.1 ug/L	0.57	101.1 ug/L	0.57	0.56%
Fe 238.204†	1504.6	10720 ug/L	45.2	10720 ug/L	45.2	0.42%
K 766.490†	19864.2	11220 ug/L	113.1	11220 ug/L	113.1	1.01%
Mg 285.213†	108214.4	16860 ug/L	145.2	16860 ug/L	145.2	0.86%
Mn 257.610†	16025.9	116.3 ug/L	1.23	116.3 ug/L	1.23	1.05%
Mo 202.031†	46.1	93.69 ug/L	5.287	93.69 ug/L	5.287	5.64%
Na 589.592†	73247.9	13420 ug/L	130.0	13420 ug/L	130.0	0.97%
Ni 231.604†	178.3	98.58 ug/L	0.938	98.58 ug/L	0.938	0.95%
Pb 220.353†	33.8	75.60 ug/L	3.802	75.60 ug/L	3.802	5.03%
Sb 206.836†	18.5	65.79 ug/L	5.826	65.79 ug/L	5.826	8.86%
Se 196.026†	54.5	641.9 ug/L	38.61	641.9 ug/L	38.61	6.01%
SiO2 251.603†	19192.4	5770 ug/L	50.3	5770 ug/L	50.3	0.87%
Sr 421.552†	1381535.2	905.0 ug/L	0.74	905.0 ug/L	0.74	0.08%
Ti 334.940†	-81.8	-0.387 ug/L	0.0501	-0.387 ug/L	0.0501	12.95%
Tl 190.801†	26.8	104.5 ug/L	11.68	104.5 ug/L	11.68	11.18%
V 290.880†	2843.0	96.39 ug/L	1.229	96.39 ug/L	1.229	1.27%
Zn 206.200†	114.3	105.4 ug/L	0.86	105.4 ug/L	0.86	0.81%

Matrix Recovery Check: 1011109-MS1

Analyte	Expected Conc.	Measured Conc.	Std. Dev.	Units	Recovery (%)
Al 396.153	10120	10810	98.911	ug/L	106.9
Ca 317.933	52570	51960	418.521	ug/L	93.9
Fe 238.204	10130	10720	45.180	ug/L	105.8
K 766.490	10750	11220	113.143	ug/L	104.6
Mg 285.213	16530	16860	145.235	ug/L	103.2
Na 589.592	12830	13420	130.007	ug/L	105.8
Ag 328.068	100.8	104.8	0.639	ug/L	104.0
As 193.696	25.59	40.03	52.990	ug/L	114.4
Ba 233.527	167.2	166.0	1.182	ug/L	98.9
Be 313.107	99.90	106.3	0.818	ug/L	106.4
Cd 214.440	100.2	105.3	0.454	ug/L	105.0
Co 228.616	100.6	101.4	1.277	ug/L	100.9
Cr 267.716	100.1	99.91	1.066	ug/L	99.8
Cu 324.752	101.1	101.1	0.565	ug/L	99.9
Mn 257.610	114.3	116.3	1.226	ug/L	102.0
Mo 202.031	107.8	93.69	5.287	ug/L	85.9
Ni 231.604	99.79	98.58	0.938	ug/L	98.8
Pb 220.353	87.99	75.60	3.802	ug/L	87.6
Sb 206.836	96.70	65.79	5.826	ug/L	69.1
Se 196.026	616.9	641.9	38.612	ug/L	105.0

SiO2 251.603	10840	5770	50.291	ug/L	-1.3
Sr 421.552	868.8	905.0	0.744	ug/L	107.2
Tl 190.801	94.62	104.5	11.680	ug/L	109.9
V 290.880	97.20	96.39	1.229	ug/L	99.2
Zn 206.200	102.7	105.4	0.859	ug/L	102.7

Sequence No.: 16
 Sample ID: 1011109-MSD1
 Analyst: Walker
 Initial Sample Wt:
 Dilution:

Autosampler Location: 32
 Date Collected: 11/29/2010 8:25:28 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: 1011109-MSD1

Analyte	Back Pressure	Flow
All	115.0 kPa	0.80 L/min

Mean Data: 1011109-MSD1

Analyte	Mean Corrected Intensity	Conc. Units	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2474381.9	101.0 %	0.34			0.34%
Sc Radial	309433.3	103.2 %	0.90			0.87%
Ag 328.068†	6207.5	103.5 ug/L	1.01	103.5 ug/L	1.01	0.97%
Al 396.153†	49650.2	10660 ug/L	63.1	10660 ug/L	63.1	0.59%
As 193.696†	1.9	30.31 ug/L	23.891	30.31 ug/L	23.891	78.81%
Ba 233.527†	1993.6	164.9 ug/L	0.47	164.9 ug/L	0.47	0.28%
Be 313.107†	28725.3	104.7 ug/L	1.02	104.7 ug/L	1.02	0.97%
B 249.677†	-118.8	-18.24 ug/L	0.609	-18.24 ug/L	0.609	3.34%
Ca 317.933†	149881.2	51980 ug/L	298.2	51980 ug/L	298.2	0.57%
Cd 214.440†	339.2	103.6 ug/L	1.34	103.6 ug/L	1.34	1.29%
Co 228.616†	420.3	99.06 ug/L	0.810	99.06 ug/L	0.810	0.82%
Cr 267.716†	616.7	98.99 ug/L	0.400	98.99 ug/L	0.400	0.40%
Cu 324.752†	9513.6	99.62 ug/L	0.995	99.62 ug/L	0.995	1.00%
Fe 238.204†	1473.5	10500 ug/L	45.7	10500 ug/L	45.7	0.44%
K 766.490†	19684.6	11120 ug/L	47.0	11120 ug/L	47.0	0.42%
Mg 285.213†	107591.1	16760 ug/L	90.8	16760 ug/L	90.8	0.54%
Mn 257.610†	15796.0	114.7 ug/L	1.28	114.7 ug/L	1.28	1.12%
Mo 202.031†	49.9	101.3 ug/L	3.86	101.3 ug/L	3.86	3.81%
Na 589.592†	72485.7	13280 ug/L	35.4	13280 ug/L	35.4	0.27%
Ni 231.604†	175.2	96.91 ug/L	1.667	96.91 ug/L	1.667	1.72%
Pb 220.353†	34.9	77.92 ug/L	4.026	77.92 ug/L	4.026	5.17%
Sb 206.836†	22.0	79.23 ug/L	5.383	79.23 ug/L	5.383	6.79%
Se 196.026†	60.3	709.9 ug/L	19.48	709.9 ug/L	19.48	2.74%
SiO2 251.603†	19308.1	5805 ug/L	44.7	5805 ug/L	44.7	0.77%
Sr 421.552†	1375706.4	901.1 ug/L	1.29	901.1 ug/L	1.29	0.14%
Ti 334.940†	-87.7	-0.415 ug/L	0.0938	-0.415 ug/L	0.0938	22.62%
Tl 190.801†	23.8	92.25 ug/L	10.420	92.25 ug/L	10.420	11.30%
V 290.880†	2789.3	94.56 ug/L	1.051	94.56 ug/L	1.051	1.11%
Zn 206.200†	112.1	103.3 ug/L	3.32	103.3 ug/L	3.32	3.21%

Matrix Recovery Check: 1011109-MSD1

Analyte	Expected Conc.	Measured Conc.	Std. Dev.	Units	Recovery (%)
Al 396.153	10120	10660	63.097	ug/L	105.4
Ca 317.933	52570	51980	298.167	ug/L	94.2
Fe 238.204	10130	10500	45.732	ug/L	103.7
K 766.490	10750	11120	47.013	ug/L	103.6
Mg 285.213	16530	16760	90.757	ug/L	102.2
Na 589.592	12830	13280	35.444	ug/L	104.4
Ag 328.068	100.8	103.5	1.005	ug/L	102.7
As 193.696	25.59	30.31	23.891	ug/L	104.7
Ba 233.527	167.2	164.9	0.466	ug/L	97.7
Be 313.107	99.90	104.7	1.015	ug/L	104.8
Cd 214.440	100.2	103.6	1.341	ug/L	103.4
Co 228.616	100.6	99.06	0.810	ug/L	98.5
Cr 267.716	100.1	98.99	0.400	ug/L	98.9
Cu 324.752	101.1	99.62	0.995	ug/L	98.5
Mn 257.610	114.3	114.7	1.282	ug/L	100.3
Mo 202.031	107.8	101.3	3.861	ug/L	93.5
Ni 231.604	99.79	96.91	1.667	ug/L	97.1
Pb 220.353	87.99	77.92	4.026	ug/L	89.9
Sb 206.836	96.70	79.23	5.383	ug/L	82.5
Se 196.026	616.9	709.9	19.482	ug/L	118.6

SiO2 251.603	10840	5805	44.660	ug/L	-0.6
Sr 421.552	868.8	901.1	1.291	ug/L	106.5
Tl 190.801	94.62	92.25	10.420	ug/L	97.6
V 290.880	97.20	94.56	1.051	ug/L	97.4
Zn 206.200	102.7	103.3	3.320	ug/L	100.7

Sequence No.: 17
 Sample ID: C101104-06
 Analyst: Walker
 Initial Sample Wt:
 Dilution:

Autosampler Location: 33
 Date Collected: 11/29/2010 8:28:34 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: C101104-06

Analyte Back Pressure Flow
 All 115.0 kPa 0.80 L/min

Mean Data: C101104-06

Analyte	Mean Corrected Intensity	Calib Conc. Units	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2515551.3	102.6 %	0.38			0.37%
Sc Radial	313647.4	104.6 %	0.60			0.57%
Ag 328.068†	7.2	0.8185 ug/L	0.33585	0.8185 ug/L	0.33585	41.03%
Al 396.153†	126.1	19.20 ug/L	4.310	19.20 ug/L	4.310	22.45%
As 193.696†	-8.7	-88.94 ug/L	22.500	-88.94 ug/L	22.500	25.30%
Ba 233.527†	813.9	67.31 ug/L	0.212	67.31 ug/L	0.212	0.31%
Be 313.107†	5.4	-0.0630 ug/L	0.02503	-0.0630 ug/L	0.02503	39.73%
B 249.677†	13.8	2.123 ug/L	0.3767	2.123 ug/L	0.3767	17.74%
Ca 317.933†	122194.5	42390 ug/L	263.8	42390 ug/L	263.8	0.62%
Cd 214.440†	-1.5	-0.4170 ug/L	1.11795	-0.4170 ug/L	1.11795	268.09%
Co 228.616†	1.2	0.4689 ug/L	0.44580	0.4689 ug/L	0.44580	95.07%
Cr 267.716†	-1.5	0.3347 ug/L	0.16912	0.3347 ug/L	0.16912	50.52%
Cu 324.752†	103.9	1.255 ug/L	0.1106	1.255 ug/L	0.1106	8.81%
Fe 238.204†	7.0	44.20 ug/L	4.265	44.20 ug/L	4.265	9.65%
K 766.490†	1131.4	631.8 ug/L	15.88	631.8 ug/L	15.88	2.51%
Mg 285.213†	41176.2	6413 ug/L	34.3	6413 ug/L	34.3	0.53%
Mn 257.610†	1977.1	14.17 ug/L	0.087	14.17 ug/L	0.087	0.62%
Mo 202.031†	5.3	10.42 ug/L	1.750	10.42 ug/L	1.750	16.80%
Na 589.592†	15060.0	2759 ug/L	10.3	2759 ug/L	10.3	0.37%
Ni 231.604†	-1.7	-0.9559 ug/L	0.81938	-0.9559 ug/L	0.81938	85.72%
Pb 220.353†	-7.6	-17.01 ug/L	5.862	-17.01 ug/L	5.862	34.47%
Sb 206.836†	1.4	3.260 ug/L	12.2665	3.260 ug/L	12.2665	376.24%
Se 196.026†	7.9	90.86 ug/L	42.350	90.86 ug/L	42.350	46.61%
SiO2 251.603†	19126.4	5768 ug/L	72.3	5768 ug/L	72.3	1.25%
Sr 421.552†	565807.1	370.6 ug/L	0.59	370.6 ug/L	0.59	0.16%
Ti 334.940†	-26.7	-0.126 ug/L	0.0141	-0.126 ug/L	0.0141	11.15%
Tl 190.801†	-1.9	-10.48 ug/L	8.764	-10.48 ug/L	8.764	83.63%
V 290.880†	-87.2	-4.428 ug/L	0.4941	-4.428 ug/L	0.4941	11.16%
Zn 206.200†	4.5	3.704 ug/L	1.0817	3.704 ug/L	1.0817	29.20%

Sequence No.: 18
 Sample ID: C101104-09
 Analyst: Walker
 Initial Sample Wt:
 Dilution:

Autosampler Location: 34
 Date Collected: 11/29/2010 8:31:39 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: C101104-09

Analyte Back Pressure Flow
 All 115.0 kPa 0.80 L/min

Mean Data: C101104-09

Analyte	Mean Corrected Intensity	Calib Conc. Units	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2493526.4	101.7 %	0.47			0.46%
Sc Radial	311541.9	103.9 %	0.30			0.29%
Ag 328.068†	25.0	1.323 ug/L	0.1770	1.323 ug/L	0.1770	13.37%
Al 396.153†	133.7	19.61 ug/L	3.685	19.61 ug/L	3.685	18.79%
As 193.696†	-8.6	-88.33 ug/L	24.208	-88.33 ug/L	24.208	27.41%
Ba 233.527†	750.8	61.83 ug/L	0.564	61.83 ug/L	0.564	0.91%
Be 313.107†	-0.2	-0.1081 ug/L	0.03847	-0.1081 ug/L	0.03847	35.61%
B 249.677†	-83.3	-12.78 ug/L	0.460	-12.78 ug/L	0.460	3.60%
Ca 317.933†	156967.3	54460 ug/L	100.4	54460 ug/L	100.4	0.18%
Cd 214.440†	-1.5	-0.3999 ug/L	0.26769	-0.3999 ug/L	0.26769	66.94%
Co 228.616†	3.1	0.9290 ug/L	0.74622	0.9290 ug/L	0.74622	80.32%
Cr 267.716†	-4.4	-0.0266 ug/L	0.19683	-0.0266 ug/L	0.19683	741.27%
Cu 324.752†	109.9	1.376 ug/L	0.1138	1.376 ug/L	0.1138	8.27%
Fe 238.204†	4.1	22.01 ug/L	17.354	22.01 ug/L	17.354	78.86%
K 766.490†	1541.5	857.8 ug/L	17.18	857.8 ug/L	17.18	2.00%
Mg 285.213†	48591.4	7567 ug/L	11.4	7567 ug/L	11.4	0.15%
Mn 257.610†	19521.4	142.2 ug/L	2.88	142.2 ug/L	2.88	2.03%
Mo 202.031†	4.5	8.706 ug/L	2.1567	8.706 ug/L	2.1567	24.77%
Na 589.592†	18419.1	3373 ug/L	11.1	3373 ug/L	11.1	0.33%
Ni 231.604†	0.2	0.1015 ug/L	1.50583	0.1015 ug/L	1.50583	>999.9%
Pb 220.353†	-6.4	-14.54 ug/L	7.311	-14.54 ug/L	7.311	50.28%
Sb 206.836†	-1.8	-9.852 ug/L	10.4279	-9.852 ug/L	10.4279	105.84%
Se 196.026†	6.0	67.62 ug/L	54.813	67.62 ug/L	54.813	81.06%
SiO2 251.603†	23879.8	7200 ug/L	83.4	7200 ug/L	83.4	1.16%
Sr 421.552†	788245.6	516.3 ug/L	1.65	516.3 ug/L	1.65	0.32%
Ti 334.940†	-35.3	-0.167 ug/L	0.0877	-0.167 ug/L	0.0877	52.54%
Tl 190.801†	-1.7	-10.63 ug/L	4.477	-10.63 ug/L	4.477	42.13%
V 290.880†	-68.8	-4.262 ug/L	0.4764	-4.262 ug/L	0.4764	11.18%
Zn 206.200†	33.0	30.45 ug/L	1.838	30.45 ug/L	1.838	6.03%

Sequence No.: 19
 Sample ID: Blank
 Analyst: Walker
 Initial Sample Wt:
 Dilution:

Autosampler Location: 35
 Date Collected: 11/29/2010 8:34:44 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: Blank

Analyte Back Pressure Flow
 All 115.0 kPa 0.80 L/min

Mean Data: Blank

Analyte	Mean Corrected Intensity	Calib Conc. Units	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2521587.3	102.9 %	0.19			0.18%
Sc Radial	307130.4	102.5 %	0.27			0.26%
Ag 328.068†	17.4	0.2860 ug/L	0.43069	0.2860 ug/L	0.43069	150.60%
Al 396.153†	4.8	1.011 ug/L	1.0625	1.011 ug/L	1.0625	105.13%
As 193.696†	-1.3	-12.80 ug/L	19.005	-12.80 ug/L	19.005	148.51%
Ba 233.527†	1.6	0.1204 ug/L	0.08896	0.1204 ug/L	0.08896	73.90%
Be 313.107†	-5.3	-0.0198 ug/L	0.05462	-0.0198 ug/L	0.05462	275.79%
B 249.677†	47.7	7.321 ug/L	0.6087	7.321 ug/L	0.6087	8.31%
Ca 317.933†	17.4	5.987 ug/L	2.8674	5.987 ug/L	2.8674	47.90%
Cd 214.440†	1.0	0.3203 ug/L	0.38324	0.3203 ug/L	0.38324	119.64%
Co 228.616†	-0.6	-0.1423 ug/L	0.93672	-0.1423 ug/L	0.93672	658.23%
Cr 267.716†	-2.0	-0.3103 ug/L	0.03504	-0.3103 ug/L	0.03504	11.29%
Cu 324.752†	69.0	0.7240 ug/L	0.25569	0.7240 ug/L	0.25569	35.32%
Fe 238.204†	2.9	20.38 ug/L	16.805	20.38 ug/L	16.805	82.44%
K 766.490†	64.0	36.82 ug/L	22.837	36.82 ug/L	22.837	62.02%
Mg 285.213†	0.6	0.0783 ug/L	1.14257	0.0783 ug/L	1.14257	>999.9%
Mn 257.610†	-11.1	-0.0832 ug/L	0.03575	-0.0832 ug/L	0.03575	42.97%
Mo 202.031†	0.1	0.1878 ug/L	1.48375	0.1878 ug/L	1.48375	790.12%
Na 589.592†	6.4	1.098 ug/L	5.9879	1.098 ug/L	5.9879	545.35%
Ni 231.604†	-0.5	-0.2716 ug/L	0.61472	-0.2716 ug/L	0.61472	226.35%
Pb 220.353†	-0.6	-1.388 ug/L	7.8196	-1.388 ug/L	7.8196	563.26%
Sb 206.836†	1.9	7.239 ug/L	11.8498	7.239 ug/L	11.8498	163.69%
Se 196.026†	1.8	21.81 ug/L	16.194	21.81 ug/L	16.194	74.24%
SiO2 251.603†	345.5	104.3 ug/L	1.58	104.3 ug/L	1.58	1.52%
Sr 421.552†	50.9	0.033 ug/L	0.0023	0.033 ug/L	0.0023	7.03%
Ti 334.940†	11.1	0.052 ug/L	0.0168	0.052 ug/L	0.0168	32.04%
Tl 190.801†	-2.0	-8.350 ug/L	2.7590	-8.350 ug/L	2.7590	33.04%
V 290.880†	-99.8	-3.537 ug/L	0.4458	-3.537 ug/L	0.4458	12.60%
Zn 206.200†	1.0	0.969 ug/L	1.0889	0.969 ug/L	1.0889	112.37%

Sequence No.: 20
 Sample ID: SEQ-CCV
 Analyst:
 Initial Sample Wt:
 Dilution:

Autosampler Location: 3
 Date Collected: 11/29/2010 8:37:49 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: SEQ-CCV

Analyte	Back Pressure	Flow
All	115.0 kPa	0.80 L/min

Mean Data: SEQ-CCV

Analyte	Mean Corrected Intensity	Calib Conc. Units	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2489303.8	101.6 %	0.75			0.74%
Sc Radial	313708.4	104.7 %	0.64			0.61%
Ag 328.068†	15651.8	258.7 ug/L	2.33	258.7 ug/L	2.33	0.90%
QC value within limits for Ag 328.068		Recovery = 103.49%				
Al 396.153†	59075.2	12680 ug/L	94.6	12680 ug/L	94.6	0.75%
QC value within limits for Al 396.153		Recovery = 101.43%				
As 193.696†	253.3	2589 ug/L	44.5	2589 ug/L	44.5	1.72%
QC value within limits for As 193.696		Recovery = 103.58%				
Ba 233.527†	6065.2	504.0 ug/L	3.36	504.0 ug/L	3.36	0.67%
QC value within limits for Ba 233.527		Recovery = 100.79%				
Be 313.107†	142759.0	520.7 ug/L	1.71	520.7 ug/L	1.71	0.33%
QC value within limits for Be 313.107		Recovery = 104.14%				
B 249.677†	33009.9	5069 ug/L	67.4	5069 ug/L	67.4	1.33%
QC value within limits for B 249.677		Recovery = 101.38%				
Ca 317.933†	36708.1	12670 ug/L	104.2	12670 ug/L	104.2	0.82%
QC value within limits for Ca 317.933		Recovery = 101.40%				
Cd 214.440†	1684.1	514.1 ug/L	3.59	514.1 ug/L	3.59	0.70%
QC value within limits for Cd 214.440		Recovery = 102.83%				
Co 228.616†	2208.1	519.7 ug/L	4.97	519.7 ug/L	4.97	0.96%
QC value within limits for Co 228.616		Recovery = 103.94%				
Cr 267.716†	15679.6	2486 ug/L	13.3	2486 ug/L	13.3	0.54%
QC value within limits for Cr 267.716		Recovery = 99.44%				
Cu 324.752†	98910.7	1032 ug/L	3.2	1032 ug/L	3.2	0.31%
QC value within limits for Cu 324.752		Recovery = 103.17%				
Fe 238.204†	1775.9	12650 ug/L	105.8	12650 ug/L	105.8	0.84%
QC value within limits for Fe 238.204		Recovery = 101.19%				
K 766.490†	44605.5	25200 ug/L	167.0	25200 ug/L	167.0	0.66%
QC value within limits for K 766.490		Recovery = 100.81%				
Mg 285.213†	82083.6	12780 ug/L	109.7	12780 ug/L	109.7	0.86%
QC value within limits for Mg 285.213		Recovery = 102.22%				
Mn 257.610†	140325.0	1023 ug/L	8.3	1023 ug/L	8.3	0.81%
QC value within limits for Mn 257.610		Recovery = 102.32%				
Mo 202.031†	253.3	513.4 ug/L	9.36	513.4 ug/L	9.36	1.82%
QC value within limits for Mo 202.031		Recovery = 102.68%				
Na 589.592†	69743.6	12700 ug/L	83.0	12700 ug/L	83.0	0.65%
QC value within limits for Na 589.592		Recovery = 101.60%				
Ni 231.604†	4577.0	2537 ug/L	13.0	2537 ug/L	13.0	0.51%
QC value within limits for Ni 231.604		Recovery = 101.48%				
Pb 220.353†	1147.2	2584 ug/L	18.3	2584 ug/L	18.3	0.71%
QC value within limits for Pb 220.353		Recovery = 103.34%				
Sb 206.836†	659.2	2507 ug/L	44.7	2507 ug/L	44.7	1.78%
QC value within limits for Sb 206.836		Recovery = 100.29%				
Se 196.026†	226.1	2668 ug/L	33.4	2668 ug/L	33.4	1.25%
QC value within limits for Se 196.026		Recovery = 106.73%				
SiO2 251.603†	33577.8	10120 ug/L	73.3	10120 ug/L	73.3	0.72%
QC value within limits for SiO2 251.603		Recovery = 101.15%				
Sr 421.552†	805617.3	527.7 ug/L	1.14	527.7 ug/L	1.14	0.22%
QC value within limits for Sr 421.552		Recovery = 105.54%				
Ti 334.940†	107363.7	507.9 ug/L	0.80	507.9 ug/L	0.80	0.16%
QC value within limits for Ti 334.940		Recovery = 101.57%				
Tl 190.801†	634.1	2593 ug/L	5.4	2593 ug/L	5.4	0.21%
QC value within limits for Tl 190.801		Recovery = 103.73%				
V 290.880†	28397.2	1001 ug/L	6.2	1001 ug/L	6.2	0.62%
QC value within limits for V 290.880		Recovery = 100.12%				

Zn 206.200†	2731.1	2575 ug/L	15.3	2575 ug/L	15.3	0.59%
-------------	--------	-----------	------	-----------	------	-------

QC value within limits for Zn 206.200 Recovery = 102.98%

All analyte(s) passed QC.

Zn 206.200†	2.0	1.916 ug/L	1.9933	1.916 ug/L	1.9933 104.03%
-------------	-----	------------	--------	------------	----------------

QC value within limits for Zn 206.200 Recovery = Not calculated
All analyte(s) passed QC.

Sequence No.: 22
 Sample ID: C101104-12
 Analyst: Walker
 Initial Sample Wt:
 Dilution:

Autosampler Location: 36
 Date Collected: 11/29/2010 8:43:59 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: C101104-12

Analyte Back Pressure Flow
 All 114.0 kPa 0.80 L/min

Mean Data: C101104-12

Analyte	Mean Corrected Intensity	Conc. Units	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2472806.4	100.9 %	0.69			0.68%
Sc Radial	309607.1	103.3 %	1.15			1.11%
Ag 328.068†	1.0	1.244 ug/L	0.2638	1.244 ug/L	0.2638	21.20%
Al 396.153†	155.5	24.62 ug/L	2.838	24.62 ug/L	2.838	11.53%
As 193.696†	-10.4	-106.5 ug/L	28.10	-106.5 ug/L	28.10	26.38%
Ba 233.527†	729.2	59.77 ug/L	0.335	59.77 ug/L	0.335	0.56%
Be 313.107†	-7.3	-0.1605 ug/L	0.02784	-0.1605 ug/L	0.02784	17.34%
B 249.677†	15.4	2.361 ug/L	0.9630	2.361 ug/L	0.9630	40.79%
Ca 317.933†	189698.8	65810 ug/L	406.2	65810 ug/L	406.2	0.62%
Cd 214.440†	1.9	0.6296 ug/L	0.41362	0.6296 ug/L	0.41362	65.69%
Co 228.616†	3.3	1.037 ug/L	1.0778	1.037 ug/L	1.0778	103.98%
Cr 267.716†	-4.3	0.1287 ug/L	0.17659	0.1287 ug/L	0.17659	137.18%
Cu 324.752†	125.4	1.628 ug/L	0.0874	1.628 ug/L	0.0874	5.37%
Fe 238.204†	4.5	23.88 ug/L	8.984	23.88 ug/L	8.984	37.63%
K 766.490†	1690.8	928.0 ug/L	23.22	928.0 ug/L	23.22	2.50%
Mg 285.213†	54385.8	8469 ug/L	13.4	8469 ug/L	13.4	0.16%
Mn 257.610†	33817.6	246.5 ug/L	2.60	246.5 ug/L	2.60	1.06%
Mo 202.031†	6.2	12.09 ug/L	3.761	12.09 ug/L	3.761	31.10%
Na 589.592†	21704.8	3973 ug/L	3.3	3973 ug/L	3.3	0.08%
Ni 231.604†	-2.7	-1.479 ug/L	0.7682	-1.479 ug/L	0.7682	51.95%
Pb 220.353†	-3.8	-8.758 ug/L	10.5667	-8.758 ug/L	10.5667	120.64%
Sb 206.836†	-1.6	-10.01 ug/L	12.592	-10.01 ug/L	12.592	125.83%
Se 196.026†	9.2	105.1 ug/L	31.79	105.1 ug/L	31.79	30.25%
SiO2 251.603†	25663.4	7735 ug/L	73.6	7735 ug/L	73.6	0.95%
Sr 421.552†	1096193.7	718.0 ug/L	0.74	718.0 ug/L	0.74	0.10%
Ti 334.940†	-28.4	-0.134 ug/L	0.0224	-0.134 ug/L	0.0224	16.73%
Tl 190.801†	-0.3	-5.957 ug/L	8.8572	-5.957 ug/L	8.8572	148.67%
V 290.880†	-68.1	-4.811 ug/L	0.2178	-4.811 ug/L	0.2178	4.53%
Zn 206.200†	187.8	176.4 ug/L	0.43	176.4 ug/L	0.43	0.25%

Sequence No.: 23
 Sample ID: C101104-15
 Analyst: Walker
 Initial Sample Wt:
 Dilution:

Autosampler Location: 37
 Date Collected: 11/29/2010 8:47:06 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: C101104-15

Analyte Back Pressure Flow
 All 117.0 kPa 0.80 L/min

Mean Data: C101104-15

Analyte	Mean Corrected Intensity	Calib Conc. Units	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2474836.5	101.0 %	0.59			0.58%
Sc Radial	311297.5	103.9 %	0.66			0.64%
Ag 328.068†	0.9	1.687 ug/L	0.2248	1.687 ug/L	0.2248	13.32%
Al 396.153†	183.0	31.35 ug/L	6.370	31.35 ug/L	6.370	20.32%
As 193.696†	-8.8	-91.12 ug/L	8.703	-91.12 ug/L	8.703	9.55%
Ba 233.527†	679.7	55.26 ug/L	0.433	55.26 ug/L	0.433	0.78%
Be 313.107†	4.5	-0.1572 ug/L	0.04158	-0.1572 ug/L	0.04158	26.45%
B 249.677†	-76.9	-11.81 ug/L	0.726	-11.81 ug/L	0.726	6.15%
Ca 317.933†	239361.1	83040 ug/L	409.2	83040 ug/L	409.2	0.49%
Cd 214.440†	4.4	1.408 ug/L	0.2999	1.408 ug/L	0.2999	21.30%
Co 228.616†	1.1	0.5579 ug/L	0.99771	0.5579 ug/L	0.99771	178.85%
Cr 267.716†	-4.8	0.2688 ug/L	0.44473	0.2688 ug/L	0.44473	165.44%
Cu 324.752†	162.9	2.139 ug/L	0.1193	2.139 ug/L	0.1193	5.58%
Fe 238.204†	3.0	10.83 ug/L	14.247	10.83 ug/L	14.247	131.57%
K 766.490†	1940.2	1048 ug/L	14.5	1048 ug/L	14.5	1.38%
Mg 285.213†	63130.7	9831 ug/L	59.3	9831 ug/L	59.3	0.60%
Mn 257.610†	54870.8	400.1 ug/L	3.16	400.1 ug/L	3.16	0.79%
Mo 202.031†	3.4	6.457 ug/L	2.3896	6.457 ug/L	2.3896	37.01%
Na 589.592†	26151.6	4785 ug/L	22.0	4785 ug/L	22.0	0.46%
Ni 231.604†	-0.3	-0.0884 ug/L	0.12612	-0.0884 ug/L	0.12612	142.72%
Pb 220.353†	-5.8	-13.63 ug/L	2.702	-13.63 ug/L	2.702	19.82%
Sb 206.836†	-1.8	-12.19 ug/L	6.102	-12.19 ug/L	6.102	50.08%
Se 196.026†	6.5	72.91 ug/L	35.467	72.91 ug/L	35.467	48.64%
SiO2 251.603†	28572.6	8607 ug/L	71.2	8607 ug/L	71.2	0.83%
Sr 421.552†	1561898.7	1023 ug/L	0.5	1023 ug/L	0.5	0.04%
Ti 334.940†	-44.4	-0.210 ug/L	0.1132	-0.210 ug/L	0.1132	53.91%
Tl 190.801†	-4.3	-23.54 ug/L	16.442	-23.54 ug/L	16.442	69.84%
V 290.880†	-55.9	-5.276 ug/L	0.4327	-5.276 ug/L	0.4327	8.20%
Zn 206.200†	414.2	389.9 ug/L	1.50	389.9 ug/L	1.50	0.38%

Sequence No.: 24
 Sample ID: C101104-18
 Analyst: Walker
 Initial Sample Wt:
 Dilution:

Autosampler Location: 38
 Date Collected: 11/29/2010 8:50:13 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: C101104-18

Analyte Back Pressure Flow
 All 116.0 kPa 0.80 L/min

Mean Data: C101104-18

Analyte	Mean Corrected Intensity	Conc. Units	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2482473.8	101.3 %	0.49			0.48%
Sc Radial	311228.6	103.8 %	1.35			1.30%
Ag 328.068†	2.2	1.265 ug/L	0.0668	1.265 ug/L	0.0668	5.28%
Al 396.153†	137.7	20.88 ug/L	5.332	20.88 ug/L	5.332	25.54%
As 193.696†	-8.7	-89.10 ug/L	5.057	-89.10 ug/L	5.057	5.68%
Ba 233.527†	728.1	59.66 ug/L	0.440	59.66 ug/L	0.440	0.74%
Be 313.107†	-0.1	-0.1349 ug/L	0.01929	-0.1349 ug/L	0.01929	14.30%
B 249.677†	-75.6	-11.60 ug/L	0.416	-11.60 ug/L	0.416	3.58%
Ca 317.933†	190192.3	65980 ug/L	196.7	65980 ug/L	196.7	0.30%
Cd 214.440†	2.5	0.8280 ug/L	0.18196	0.8280 ug/L	0.18196	21.98%
Co 228.616†	4.0	1.196 ug/L	0.5761	1.196 ug/L	0.5761	48.18%
Cr 267.716†	-3.5	0.2651 ug/L	0.31402	0.2651 ug/L	0.31402	118.46%
Cu 324.752†	143.4	1.818 ug/L	0.0498	1.818 ug/L	0.0498	2.74%
Fe 238.204†	2.1	6.081 ug/L	11.3315	6.081 ug/L	11.3315	186.35%
K 766.490†	1739.0	955.4 ug/L	25.35	955.4 ug/L	25.35	2.65%
Mg 285.213†	54552.3	8495 ug/L	20.5	8495 ug/L	20.5	0.24%
Mn 257.610†	33571.4	244.7 ug/L	1.11	244.7 ug/L	1.11	0.45%
Mo 202.031†	5.3	10.21 ug/L	3.964	10.21 ug/L	3.964	38.82%
Na 589.592†	21585.5	3951 ug/L	9.5	3951 ug/L	9.5	0.24%
Ni 231.604†	1.1	0.5990 ug/L	0.23562	0.5990 ug/L	0.23562	39.34%
Pb 220.353†	-7.5	-17.23 ug/L	8.652	-17.23 ug/L	8.652	50.20%
Sb 206.836†	2.6	6.076 ug/L	3.1496	6.076 ug/L	3.1496	51.83%
Se 196.026†	12.1	139.5 ug/L	26.16	139.5 ug/L	26.16	18.75%
SiO2 251.603†	25493.4	7684 ug/L	44.8	7684 ug/L	44.8	0.58%
Sr 421.552†	1095886.0	717.8 ug/L	1.66	717.8 ug/L	1.66	0.23%
Ti 334.940†	-37.5	-0.178 ug/L	0.0712	-0.178 ug/L	0.0712	40.12%
Tl 190.801†	-3.5	-18.83 ug/L	11.949	-18.83 ug/L	11.949	63.45%
V 290.880†	-86.6	-5.489 ug/L	0.3151	-5.489 ug/L	0.3151	5.74%
Zn 206.200†	184.1	172.9 ug/L	0.89	172.9 ug/L	0.89	0.52%

Sequence No.: 25
 Sample ID: C101104-21
 Analyst: Walker
 Initial Sample Wt:
 Dilution:

Autosampler Location: 39
 Date Collected: 11/29/2010 8:53:18 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: C101104-21

Analyte Back Pressure Flow
 All 115.0 kPa 0.80 L/min

Mean Data: C101104-21

Analyte	Mean Corrected Intensity	Conc. Units	Calib	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2470282.1	100.8 %		0.39			0.38%
Sc Radial	308556.3	102.9 %		0.47			0.46%
Ag 328.068†	4.0	0.8245 ug/L		0.32496	0.8245 ug/L	0.32496	39.41%
Al 396.153†	169.7	20.60 ug/L		6.211	20.60 ug/L	6.211	30.15%
As 193.696†	-10.4	-106.5 ug/L		29.86	-106.5 ug/L	29.86	28.04%
Ba 233.527†	628.0	51.60 ug/L		0.287	51.60 ug/L	0.287	0.56%
Be 313.107†	-1.8	-0.1239 ug/L		0.01838	-0.1239 ug/L	0.01838	14.83%
B 249.677†	-144.6	-22.20 ug/L		0.438	-22.20 ug/L	0.438	1.97%
Ca 317.933†	196390.1	68140 ug/L		162.2	68140 ug/L	162.2	0.24%
Cd 214.440†	0.5	0.2122 ug/L		0.16702	0.2122 ug/L	0.16702	78.69%
Co 228.616†	5.3	1.475 ug/L		0.4289	1.475 ug/L	0.4289	29.07%
Cr 267.716†	-7.3	-0.3990 ug/L		0.39627	-0.3990 ug/L	0.39627	99.32%
Cu 324.752†	122.4	1.387 ug/L		0.0968	1.387 ug/L	0.0968	6.98%
Fe 238.204†	14.1	93.71 ug/L		20.028	93.71 ug/L	20.028	21.37%
K 766.490†	1462.4	824.7 ug/L		9.52	824.7 ug/L	9.52	1.15%
Mg 285.213†	50879.7	7924 ug/L		24.2	7924 ug/L	24.2	0.30%
Mn 257.610†	5043.4	36.51 ug/L		0.284	36.51 ug/L	0.284	0.78%
Mo 202.031†	3.7	6.582 ug/L		1.1700	6.582 ug/L	1.1700	17.78%
Na 589.592†	12781.8	2339 ug/L		16.6	2339 ug/L	16.6	0.71%
Ni 231.604†	-1.5	-1.049 ug/L		1.9179	-1.049 ug/L	1.9179	182.85%
Pb 220.353†	-9.8	-20.97 ug/L		10.236	-20.97 ug/L	10.236	48.81%
Sb 206.836†	-0.9	-5.794 ug/L		10.1510	-5.794 ug/L	10.1510	175.20%
Se 196.026†	10.5	119.8 ug/L		40.15	119.8 ug/L	40.15	33.52%
SiO2 251.603†	22360.9	6747 ug/L		32.3	6747 ug/L	32.3	0.48%
Sr 421.552†	546004.5	357.7 ug/L		0.48	357.7 ug/L	0.48	0.13%
Ti 334.940†	-42.9	-0.203 ug/L		0.0435	-0.203 ug/L	0.0435	21.47%
Tl 190.801†	-1.4	-8.910 ug/L		14.9418	-8.910 ug/L	14.9418	167.70%
V 290.880†	-62.9	-3.956 ug/L		0.3781	-3.956 ug/L	0.3781	9.56%
Zn 206.200†	18.4	17.12 ug/L		2.340	17.12 ug/L	2.340	13.67%

Sequence No.: 26
 Sample ID: C101104-24
 Analyst: Walker
 Initial Sample Wt:
 Dilution:

Autosampler Location: 40
 Date Collected: 11/29/2010 8:56:22 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: C101104-24

Analyte Back Pressure Flow
 All 117.0 kPa 0.80 L/min

Mean Data: C101104-24

Analyte	Mean Corrected Intensity	Conc. Units	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2482681.6	101.3 %	0.32			0.32%
Sc Radial	310643.6	103.6 %	0.70			0.68%
Ag 328.068†	31.6	1.380 ug/L	0.2660	1.380 ug/L	0.2660	19.27%
Al 396.153†	163.0	23.52 ug/L	4.174	23.52 ug/L	4.174	17.75%
As 193.696†	-9.4	-94.38 ug/L	20.585	-94.38 ug/L	20.585	21.81%
Ba 233.527†	588.8	47.87 ug/L	0.486	47.87 ug/L	0.486	1.01%
Be 313.107†	-11.5	-0.1956 ug/L	0.00942	-0.1956 ug/L	0.00942	4.82%
B 249.677†	-170.4	-26.16 ug/L	0.458	-26.16 ug/L	0.458	1.75%
Ca 317.933†	212049.1	73570 ug/L	504.7	73570 ug/L	504.7	0.69%
Cd 214.440†	-0.7	-0.1724 ug/L	0.32261	-0.1724 ug/L	0.32261	187.14%
Co 228.616†	9.0	2.256 ug/L	0.3700	2.256 ug/L	0.3700	16.40%
Cr 267.716†	-3.5	-1.384 ug/L	0.5108	-1.384 ug/L	0.5108	36.91%
Cu 324.752†	122.5	1.627 ug/L	0.3425	1.627 ug/L	0.3425	21.05%
Fe 238.204†	745.3	5311 ug/L	19.6	5311 ug/L	19.6	0.37%
K 766.490†	2172.9	1204 ug/L	27.0	1204 ug/L	27.0	2.24%
Mg 285.213†	47000.0	7318 ug/L	28.7	7318 ug/L	28.7	0.39%
Mn 257.610†	428834.0	3130 ug/L	3.5	3130 ug/L	3.5	0.11%
Mo 202.031†	4.8	9.449 ug/L	2.6444	9.449 ug/L	2.6444	27.99%
Na 589.592†	15649.2	2861 ug/L	17.2	2861 ug/L	17.2	0.60%
Ni 231.604†	-0.9	-0.8169 ug/L	0.39781	-0.8169 ug/L	0.39781	48.70%
Pb 220.353†	-10.0	-22.72 ug/L	9.442	-22.72 ug/L	9.442	41.55%
Sb 206.836†	-1.1	-8.442 ug/L	10.4804	-8.442 ug/L	10.4804	124.15%
Se 196.026†	11.6	132.1 ug/L	66.96	132.1 ug/L	66.96	50.70%
SiO2 251.603†	31388.3	9460 ug/L	53.9	9460 ug/L	53.9	0.57%
Sr 421.552†	1085310.6	710.9 ug/L	0.97	710.9 ug/L	0.97	0.14%
Ti 334.940†	-41.4	-0.196 ug/L	0.0100	-0.196 ug/L	0.0100	5.10%
Tl 190.801†	-5.4	-31.19 ug/L	13.512	-31.19 ug/L	13.512	43.32%
V 290.880†	-50.1	-5.136 ug/L	0.2440	-5.136 ug/L	0.2440	4.75%
Zn 206.200†	6.9	5.089 ug/L	2.0605	5.089 ug/L	2.0605	40.49%

Sequence No.: 27
 Sample ID: C101104-27
 Analyst: Walker
 Initial Sample Wt:
 Dilution:

Autosampler Location: 41
 Date Collected: 11/29/2010 8:59:28 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: C101104-27

Analyte Back Pressure Flow
 All 114.0 kPa 0.80 L/min

Mean Data: C101104-27

Analyte	Mean Corrected Intensity	Calib Conc. Units	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2474972.0	101.0 %	0.83			0.82%
Sc Radial	309229.2	103.2 %	0.50			0.49%
Ag 328.068†	-8.0	0.6683 ug/L	0.19541	0.6683 ug/L	0.19541	29.24%
Al 396.153†	129.2	14.29 ug/L	0.828	14.29 ug/L	0.828	5.80%
As 193.696†	-9.2	-94.61 ug/L	17.847	-94.61 ug/L	17.847	18.86%
Ba 233.527†	549.3	45.04 ug/L	0.280	45.04 ug/L	0.280	0.62%
Be 313.107†	-7.3	-0.1371 ug/L	0.04422	-0.1371 ug/L	0.04422	32.25%
B 249.677†	-148.3	-22.77 ug/L	0.832	-22.77 ug/L	0.832	3.66%
Ca 317.933†	178376.8	61890 ug/L	304.9	61890 ug/L	304.9	0.49%
Cd 214.440†	1.4	0.4618 ug/L	0.08633	0.4618 ug/L	0.08633	18.70%
Co 228.616†	16.9	4.178 ug/L	0.9650	4.178 ug/L	0.9650	23.10%
Cr 267.716†	-2.0	0.4194 ug/L	0.66540	0.4194 ug/L	0.66540	158.64%
Cu 324.752†	209.8	2.325 ug/L	0.1818	2.325 ug/L	0.1818	7.82%
Fe 238.204†	4.7	27.31 ug/L	16.085	27.31 ug/L	16.085	58.89%
K 766.490†	1386.7	776.4 ug/L	26.53	776.4 ug/L	26.53	3.42%
Mg 285.213†	48427.6	7542 ug/L	8.9	7542 ug/L	8.9	0.12%
Mn 257.610†	1417.1	10.04 ug/L	0.121	10.04 ug/L	0.121	1.20%
Mo 202.031†	8.0	15.44 ug/L	1.578	15.44 ug/L	1.578	10.22%
Na 589.592†	12868.9	2354 ug/L	21.3	2354 ug/L	21.3	0.90%
Ni 231.604†	7.2	3.839 ug/L	1.6795	3.839 ug/L	1.6795	43.75%
Pb 220.353†	-7.7	-16.62 ug/L	8.658	-16.62 ug/L	8.658	52.10%
Sb 206.836†	0.9	0.9924 ug/L	1.58906	0.9924 ug/L	1.58906	160.12%
Se 196.026†	10.1	115.0 ug/L	81.82	115.0 ug/L	81.82	71.14%
SiO2 251.603†	20251.1	6109 ug/L	72.8	6109 ug/L	72.8	1.19%
Sr 421.552†	602004.6	394.3 ug/L	0.69	394.3 ug/L	0.69	0.18%
Ti 334.940†	-19.8	-0.094 ug/L	0.0320	-0.094 ug/L	0.0320	34.17%
Tl 190.801†	-2.1	-11.82 ug/L	20.787	-11.82 ug/L	20.787	175.80%
V 290.880†	-88.3	-4.767 ug/L	0.8790	-4.767 ug/L	0.8790	18.44%
Zn 206.200†	139.2	131.1 ug/L	1.70	131.1 ug/L	1.70	1.30%

Sequence No.: 28
 Sample ID: C101104-30
 Analyst: Walker
 Initial Sample Wt:
 Dilution:

Autosampler Location: 42
 Date Collected: 11/29/2010 9:02:32 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: C101104-30

Analyte Back Pressure Flow
 All 115.0 kPa 0.80 L/min

Mean Data: C101104-30

Analyte	Mean Corrected Intensity	Conc. Units	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2416877.2	98.61 %	0.450			0.46%
Sc Radial	307521.5	102.6 %	0.39			0.38%
Ag 328.068†	-12.3	0.3725 ug/L	0.12777	0.3725 ug/L	0.12777	34.30%
Al 396.153†	444.3	26.76 ug/L	0.725	26.76 ug/L	0.725	2.71%
As 193.696†	-9.2	-96.19 ug/L	13.877	-96.19 ug/L	13.877	14.43%
Ba 233.527†	279.7	21.44 ug/L	0.222	21.44 ug/L	0.222	1.04%
Be 313.107†	-11.2	-0.3649 ug/L	0.02864	-0.3649 ug/L	0.02864	7.85%
B 249.677†	-286.7	-44.03 ug/L	0.741	-44.03 ug/L	0.741	1.68%
Ca 317.933†	646540.1	224300 ug/L	961.1	224300 ug/L	961.1	0.43%
Cd 214.440†	-4.3	-1.200 ug/L	0.2598	-1.200 ug/L	0.2598	21.65%
Co 228.616†	6.4	2.098 ug/L	0.9720	2.098 ug/L	0.9720	46.34%
Cr 267.716†	-6.1	0.2356 ug/L	0.45916	0.2356 ug/L	0.45916	194.89%
Cu 324.752†	124.0	0.8738 ug/L	0.20384	0.8738 ug/L	0.20384	23.33%
Fe 238.204†	332.5	2359 ug/L	11.4	2359 ug/L	11.4	0.48%
K 766.490†	3522.4	2032 ug/L	32.1	2032 ug/L	32.1	1.58%
Mg 285.213†	129001.2	20090 ug/L	65.3	20090 ug/L	65.3	0.33%
Mn 257.610†	159013.4	1160 ug/L	0.9	1160 ug/L	0.9	0.07%
Mo 202.031†	5.2	6.890 ug/L	2.2561	6.890 ug/L	2.2561	32.74%
Na 589.592†	53438.8	9803 ug/L	38.8	9803 ug/L	38.8	0.40%
Ni 231.604†	0.2	-1.576 ug/L	1.9146	-1.576 ug/L	1.9146	121.49%
Pb 220.353†	-12.6	-20.37 ug/L	6.967	-20.37 ug/L	6.967	34.20%
Sb 206.836†	-1.8	-7.831 ug/L	10.1724	-7.831 ug/L	10.1724	129.90%
Se 196.026†	14.6	154.0 ug/L	90.15	154.0 ug/L	90.15	58.55%
SiO2 251.603†	37331.7	11290 ug/L	41.0	11290 ug/L	41.0	0.36%
Sr 421.552†	Saturated2					
Ti 334.940†	-53.7	-0.254 ug/L	0.0523	-0.254 ug/L	0.0523	20.61%
Tl 190.801†	-1.7	-15.06 ug/L	14.318	-15.06 ug/L	14.318	95.07%
V 290.880†	-15.2	-4.398 ug/L	0.2401	-4.398 ug/L	0.2401	5.46%
Zn 206.200†	9.5	10.93 ug/L	0.968	10.93 ug/L	0.968	8.86%

Sequence No.: 29
 Sample ID: C101104-33
 Analyst: Walker
 Initial Sample Wt:
 Dilution:

Autosampler Location: 43
 Date Collected: 11/29/2010 9:06:37 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: C101104-33

Analyte Back Pressure Flow
 All 115.0 kPa 0.80 L/min

Mean Data: C101104-33

Analyte	Mean Corrected Intensity	Conc. Units	Calib	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2392879.0	97.64 %		0.440			0.45%
Sc Radial	303592.1	101.3 %		0.59			0.58%
Ag 328.068†	-24.0	0.3347 ug/L		0.13506	0.3347 ug/L	0.13506	40.36%
Al 396.153†	447.7	24.15 ug/L		5.206	24.15 ug/L	5.206	21.56%
As 193.696†	-8.1	-86.89 ug/L		24.233	-86.89 ug/L	24.233	27.89%
Ba 233.527†	233.0	17.45 ug/L		0.188	17.45 ug/L	0.188	1.08%
Be 313.107†	-10.9	-0.3796 ug/L		0.01887	-0.3796 ug/L	0.01887	4.97%
B 249.677†	-313.0	-48.07 ug/L		0.786	-48.07 ug/L	0.786	1.63%
Ca 317.933†	686181.3	238100 ug/L		1232.7	238100 ug/L	1232.7	0.52%
Cd 214.440†	8.2	2.634 ug/L		0.6778	2.634 ug/L	0.6778	25.73%
Co 228.616†	2.7	1.324 ug/L		0.9083	1.324 ug/L	0.9083	68.61%
Cr 267.716†	-6.3	0.8965 ug/L		0.33516	0.8965 ug/L	0.33516	37.39%
Cu 324.752†	167.5	1.298 ug/L		0.3066	1.298 ug/L	0.3066	23.62%
Fe 238.204†	18.2	115.6 ug/L		8.93	115.6 ug/L	8.93	7.72%
K 766.490†	2938.7	1700 ug/L		16.2	1700 ug/L	16.2	0.96%
Mg 285.213†	142707.6	22220 ug/L		163.4	22220 ug/L	163.4	0.74%
Mn 257.610†	27407.8	199.4 ug/L		1.01	199.4 ug/L	1.01	0.51%
Mo 202.031†	2.2	0.4082 ug/L		5.05952	0.4082 ug/L	5.05952	>999.9%
Na 589.592†	58780.9	10780 ug/L		65.5	10780 ug/L	65.5	0.61%
Ni 231.604†	1.3	-0.9840 ug/L		1.04959	-0.9840 ug/L	1.04959	106.67%
Pb 220.353†	-8.3	-10.16 ug/L		6.040	-10.16 ug/L	6.040	59.46%
Sb 206.836†	-0.3	-2.501 ug/L		12.7004	-2.501 ug/L	12.7004	507.83%
Se 196.026†	19.5	211.6 ug/L		58.66	211.6 ug/L	58.66	27.72%
SiO2 251.603†	36085.6	10920 ug/L		98.4	10920 ug/L	98.4	0.90%
Sr 421.552†	Saturated2						
Ti 334.940†	-70.1	-0.332 ug/L		0.0909	-0.332 ug/L	0.0909	27.38%
Tl 190.801†	-0.4	-8.402 ug/L		13.3368	-8.402 ug/L	13.3368	158.73%
V 290.880†	-4.3	-3.977 ug/L		0.6640	-3.977 ug/L	0.6640	16.69%
Zn 206.200†	611.7	580.3 ug/L		6.54	580.3 ug/L	6.54	1.13%

Sequence No.: 30
 Sample ID: C101104-36
 Analyst: Walker
 Initial Sample Wt:
 Dilution:

Autosampler Location: 44
 Date Collected: 11/29/2010 9:10:48 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: C101104-36

Analyte Back Pressure Flow
 All 115.0 kPa 0.80 L/min

Mean Data: C101104-36

Analyte	Mean Corrected Intensity	Conc. Units	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2466457.0	100.6 %	0.91			0.90%
Sc Radial	306512.6	102.3 %	0.23			0.23%
Ag 328.068†	-2.0	0.7358 ug/L	0.48957	0.7358 ug/L	0.48957	66.54%
Al 396.153†	132.8	13.09 ug/L	2.729	13.09 ug/L	2.729	20.84%
As 193.696†	-10.4	-106.3 ug/L	41.45	-106.3 ug/L	41.45	38.98%
Ba 233.527†	683.2	56.21 ug/L	0.506	56.21 ug/L	0.506	0.90%
Be 313.107†	-8.7	-0.1481 ug/L	0.05141	-0.1481 ug/L	0.05141	34.72%
B 249.677†	-162.2	-24.90 ug/L	0.363	-24.90 ug/L	0.363	1.46%
Ca 317.933†	193720.0	67210 ug/L	68.0	67210 ug/L	68.0	0.10%
Cd 214.440†	0.3	0.1306 ug/L	0.14322	0.1306 ug/L	0.14322	109.67%
Co 228.616†	4.2	1.219 ug/L	0.6480	1.219 ug/L	0.6480	53.14%
Cr 267.716†	-5.1	-0.0459 ug/L	0.55167	-0.0459 ug/L	0.55167	>999.9%
Cu 324.752†	128.7	1.457 ug/L	0.0814	1.457 ug/L	0.0814	5.59%
Fe 238.204†	5.6	33.18 ug/L	9.896	33.18 ug/L	9.896	29.83%
K 766.490†	1408.1	793.3 ug/L	33.50	793.3 ug/L	33.50	4.22%
Mg 285.213†	51212.5	7975 ug/L	19.0	7975 ug/L	19.0	0.24%
Mn 257.610†	1396.8	9.884 ug/L	0.1009	9.884 ug/L	0.1009	1.02%
Mo 202.031†	3.8	6.808 ug/L	1.0396	6.808 ug/L	1.0396	15.27%
Na 589.592†	12882.6	2357 ug/L	8.8	2357 ug/L	8.8	0.37%
Ni 231.604†	0.1	-0.1393 ug/L	0.84890	-0.1393 ug/L	0.84890	609.38%
Pb 220.353†	-8.3	-17.78 ug/L	13.912	-17.78 ug/L	13.912	78.26%
Sb 206.836†	-0.7	-5.291 ug/L	2.8974	-5.291 ug/L	2.8974	54.76%
Se 196.026†	8.8	100.1 ug/L	13.60	100.1 ug/L	13.60	13.59%
SiO2 251.603†	22564.6	6808 ug/L	69.4	6808 ug/L	69.4	1.02%
Sr 421.552†	561183.5	367.6 ug/L	0.20	367.6 ug/L	0.20	0.05%
Ti 334.940†	-42.6	-0.202 ug/L	0.0305	-0.202 ug/L	0.0305	15.11%
Tl 190.801†	-2.8	-14.73 ug/L	9.395	-14.73 ug/L	9.395	63.80%
V 290.880†	-92.3	-4.980 ug/L	0.7071	-4.980 ug/L	0.7071	14.20%
Zn 206.200†	7.1	6.430 ug/L	1.9939	6.430 ug/L	1.9939	31.01%

Sequence No.: 31
 Sample ID: Blank
 Analyst: Walker
 Initial Sample Wt:
 Dilution:

Autosampler Location: 45
 Date Collected: 11/29/2010 9:13:53 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: Blank

Analyte	Back Pressure	Flow
All	115.0 kPa	0.80 L/min

Mean Data: Blank

Analyte	Mean Corrected Intensity	Conc. Units	Calib	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2510515.3	102.4 %		0.57			0.55%
Sc Radial	306409.2	102.2 %		0.79			0.77%
Ag 328.068†	21.4	0.3490 ug/L		0.02228	0.3490 ug/L	0.02228	6.38%
Al 396.153†	43.6	9.228 ug/L		2.9456	9.228 ug/L	2.9456	31.92%
As 193.696†	0.1	0.8155 ug/L		20.74079	0.8155 ug/L	20.74079	>999.9%
Ba 233.527†	3.8	0.3057 ug/L		0.15660	0.3057 ug/L	0.15660	51.23%
Be 313.107†	4.7	0.0165 ug/L		0.00794	0.0165 ug/L	0.00794	48.01%
B 249.677†	122.9	18.88 ug/L		0.253	18.88 ug/L	0.253	1.34%
Ca 317.933†	34.6	11.97 ug/L		2.424	11.97 ug/L	2.424	20.26%
Cd 214.440†	0.6	0.1864 ug/L		0.23949	0.1864 ug/L	0.23949	128.45%
Co 228.616†	0.5	0.1117 ug/L		0.21449	0.1117 ug/L	0.21449	191.96%
Cr 267.716†	-0.9	-0.1485 ug/L		0.19928	-0.1485 ug/L	0.19928	134.21%
Cu 324.752†	88.8	0.9292 ug/L		0.36043	0.9292 ug/L	0.36043	38.79%
Fe 238.204†	1.9	13.37 ug/L		5.067	13.37 ug/L	5.067	37.89%
K 766.490†	64.6	38.18 ug/L		15.383	38.18 ug/L	15.383	40.29%
Mg 285.213†	4.5	0.6886 ug/L		0.30201	0.6886 ug/L	0.30201	43.86%
Mn 257.610†	-12.9	-0.0960 ug/L		0.06615	-0.0960 ug/L	0.06615	68.91%
Mo 202.031†	-0.6	-1.191 ug/L		0.7883	-1.191 ug/L	0.7883	66.20%
Na 589.592†	103.5	19.05 ug/L		1.192	19.05 ug/L	1.192	6.26%
Ni 231.604†	1.1	0.6068 ug/L		0.32760	0.6068 ug/L	0.32760	53.99%
Pb 220.353†	-2.6	-5.924 ug/L		5.2669	-5.924 ug/L	5.2669	88.91%
Sb 206.836†	-1.3	-4.913 ug/L		17.9038	-4.913 ug/L	17.9038	364.40%
Se 196.026†	0.6	6.632 ug/L		28.1217	6.632 ug/L	28.1217	424.00%
SiO2 251.603†	1189.7	359.0 ug/L		2.19	359.0 ug/L	2.19	0.61%
Sr 421.552†	80.8	0.053 ug/L		0.0092	0.053 ug/L	0.0092	17.41%
Ti 334.940†	19.4	0.092 ug/L		0.0038	0.092 ug/L	0.0038	4.14%
Tl 190.801†	-0.5	-2.076 ug/L		11.3368	-2.076 ug/L	11.3368	546.09%
V 290.880†	-127.7	-4.526 ug/L		0.2137	-4.526 ug/L	0.2137	4.72%
Zn 206.200†	0.7	0.622 ug/L		0.3655	0.622 ug/L	0.3655	58.75%

Sequence No.: 32
 Sample ID: SEQ-CCV
 Analyst:
 Initial Sample Wt:
 Dilution:

Autosampler Location: 3
 Date Collected: 11/29/2010 9:16:57 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: SEQ-CCV

Analyte	Back Pressure	Flow
All	115.0 kPa	0.80 L/min

Mean Data: SEQ-CCV

Analyte	Mean Corrected Intensity	Calib Conc. Units	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2484150.1	101.4 %	0.55			0.54%
Sc Radial	313395.7	104.6 %	0.90			0.86%
Ag 328.068†	15711.8	259.7 ug/L	3.23	259.7 ug/L	3.23	1.24%
QC value within limits for Ag 328.068		Recovery = 103.87%				
Al 396.153†	59315.6	12730 ug/L	65.2	12730 ug/L	65.2	0.51%
QC value within limits for Al 396.153		Recovery = 101.84%				
As 193.696†	247.1	2527 ug/L	16.3	2527 ug/L	16.3	0.64%
QC value within limits for As 193.696		Recovery = 101.06%				
Ba 233.527†	6120.5	508.6 ug/L	7.12	508.6 ug/L	7.12	1.40%
QC value within limits for Ba 233.527		Recovery = 101.72%				
Be 313.107†	143355.9	522.9 ug/L	1.83	522.9 ug/L	1.83	0.35%
QC value within limits for Be 313.107		Recovery = 104.58%				
B 249.677†	33156.2	5091 ug/L	65.1	5091 ug/L	65.1	1.28%
QC value within limits for B 249.677		Recovery = 101.83%				
Ca 317.933†	36948.4	12760 ug/L	44.5	12760 ug/L	44.5	0.35%
QC value within limits for Ca 317.933		Recovery = 102.07%				
Cd 214.440†	1674.2	511.1 ug/L	2.08	511.1 ug/L	2.08	0.41%
QC value within limits for Cd 214.440		Recovery = 102.22%				
Co 228.616†	2192.0	515.9 ug/L	2.85	515.9 ug/L	2.85	0.55%
QC value within limits for Co 228.616		Recovery = 103.18%				
Cr 267.716†	15618.2	2476 ug/L	36.0	2476 ug/L	36.0	1.45%
QC value within limits for Cr 267.716		Recovery = 99.05%				
Cu 324.752†	99278.4	1036 ug/L	1.3	1036 ug/L	1.3	0.13%
QC value within limits for Cu 324.752		Recovery = 103.56%				
Fe 238.204†	1765.4	12570 ug/L	100.0	12570 ug/L	100.0	0.80%
QC value within limits for Fe 238.204		Recovery = 100.59%				
K 766.490†	44675.4	25240 ug/L	79.1	25240 ug/L	79.1	0.31%
QC value within limits for K 766.490		Recovery = 100.97%				
Mg 285.213†	82283.6	12810 ug/L	58.7	12810 ug/L	58.7	0.46%
QC value within limits for Mg 285.213		Recovery = 102.47%				
Mn 257.610†	140911.1	1028 ug/L	12.5	1028 ug/L	12.5	1.22%
QC value within limits for Mn 257.610		Recovery = 102.75%				
Mo 202.031†	252.3	511.4 ug/L	9.39	511.4 ug/L	9.39	1.84%
QC value within limits for Mo 202.031		Recovery = 102.29%				
Na 589.592†	69981.1	12740 ug/L	59.5	12740 ug/L	59.5	0.47%
QC value within limits for Na 589.592		Recovery = 101.95%				
Ni 231.604†	4528.7	2510 ug/L	5.4	2510 ug/L	5.4	0.22%
QC value within limits for Ni 231.604		Recovery = 100.40%				
Pb 220.353†	1143.8	2576 ug/L	18.6	2576 ug/L	18.6	0.72%
QC value within limits for Pb 220.353		Recovery = 103.04%				
Sb 206.836†	648.5	2466 ug/L	52.3	2466 ug/L	52.3	2.12%
QC value within limits for Sb 206.836		Recovery = 98.65%				
Se 196.026†	220.0	2597 ug/L	6.8	2597 ug/L	6.8	0.26%
QC value within limits for Se 196.026		Recovery = 103.88%				
SiO2 251.603†	33730.3	10160 ug/L	106.0	10160 ug/L	106.0	1.04%
QC value within limits for SiO2 251.603		Recovery = 101.61%				
Sr 421.552†	812897.8	532.5 ug/L	0.95	532.5 ug/L	0.95	0.18%
QC value within limits for Sr 421.552		Recovery = 106.50%				
Ti 334.940†	107328.5	507.7 ug/L	1.43	507.7 ug/L	1.43	0.28%
QC value within limits for Ti 334.940		Recovery = 101.54%				
Tl 190.801†	635.3	2598 ug/L	16.4	2598 ug/L	16.4	0.63%
QC value within limits for Tl 190.801		Recovery = 103.91%				
V 290.880†	28365.7	1000 ug/L	14.3	1000 ug/L	14.3	1.43%
QC value within limits for V 290.880		Recovery = 100.01%				

Zn 206.200†	2718.1	2562 ug/L	16.6	2562 ug/L	16.6	0.65%
-------------	--------	-----------	------	-----------	------	-------

QC value within limits for Zn 206.200 Recovery = 102.49%

All analyte(s) passed QC.

Sequence No.: 33
Sample ID: SEQ-CCB
Analyst:
Initial Sample Wt:
Dilution:

Autosampler Location: 1.
Date Collected: 11/29/2010 9:20:03 AM
Data Type: Original
Initial Sample Vol:
Sample Prep Vol:

Nebulizer Parameters: SEO-CCB

Analyte	Back Pressure	Flow
All	115.0 kPa	0.80 L/min

Mean Data: SEQ-CCB

		Mean Corrected		Calib		Sample			
Analyte		Intensity	Conc.	Units	Std.Dev.	Conc.	Units	Std.Dev.	RSD
Sc Axial		2482801.4	101.3	%	0.32				0.32%
Sc Radial		307997.4	102.8	%	0.28				0.27%
Ag	328.068†	22.1	0.3642	ug/L	0.25975	0.3642	ug/L	0.25975	71.32%
	QC value within limits for Ag	328.068	Recovery =	Not calculated					
Al	396.153†	3.5	0.5598	ug/L	3.40445	0.5598	ug/L	3.40445	608.10%
	QC value within limits for Al	396.153	Recovery =	Not calculated					
As	193.696†	-1.2	-11.92	ug/L	17.155	-11.92	ug/L	17.155	143.90%
	QC value within limits for As	193.696	Recovery =	Not calculated					
Ba	233.527†	2.4	0.1757	ug/L	0.09086	0.1757	ug/L	0.09086	51.73%
	QC value within limits for Ba	233.527	Recovery =	Not calculated					
Be	313.107†	5.9	0.0203	ug/L	0.07846	0.0203	ug/L	0.07846	387.42%
	QC value within limits for Be	313.107	Recovery =	Not calculated					
B	249.677†	289.3	44.42	ug/L	3.903	44.42	ug/L	3.903	8.79%
	QC value within limits for B	249.677	Recovery =	Not calculated					
Ca	317.933†	6.5	2.123	ug/L	0.4201	2.123	ug/L	0.4201	19.79%
	QC value within limits for Ca	317.933	Recovery =	Not calculated					
Cd	214.440†	1.7	0.5261	ug/L	0.37076	0.5261	ug/L	0.37076	70.48%
	QC value within limits for Cd	214.440	Recovery =	Not calculated					
Co	228.616†	-3.2	-0.7596	ug/L	0.92721	-0.7596	ug/L	0.92721	122.06%
	QC value within limits for Co	228.616	Recovery =	Not calculated					
Cr	267.716†	1.2	0.1902	ug/L	0.55127	0.1902	ug/L	0.55127	289.89%
	QC value within limits for Cr	267.716	Recovery =	Not calculated					
Cu	324.752†	115.3	1.206	ug/L	0.1170	1.206	ug/L	0.1170	9.69%
	QC value within limits for Cu	324.752	Recovery =	Not calculated					
Fe	238.204†	2.4	17.42	ug/L	8.542	17.42	ug/L	8.542	49.04%
	QC value within limits for Fe	238.204	Recovery =	Not calculated					
K	766.490†	72.4	40.26	ug/L	14.072	40.26	ug/L	14.072	34.95%
	QC value within limits for K	766.490	Recovery =	Not calculated					
Mg	285.213†	8.7	1.305	ug/L	1.0189	1.305	ug/L	1.0189	78.08%
	QC value within limits for Mg	285.213	Recovery =	Not calculated					
Mn	257.610†	15.1	0.1053	ug/L	0.04477	0.1053	ug/L	0.04477	42.51%
	QC value within limits for Mn	257.610	Recovery =	Not calculated					
Mo	202.031†	2.7	5.527	ug/L	1.0070	5.527	ug/L	1.0070	18.22%
	QC value within limits for Mo	202.031	Recovery =	Not calculated					
Na	589.592†	11.0	1.783	ug/L	2.6551	1.783	ug/L	2.6551	148.94%
	QC value within limits for Na	589.592	Recovery =	Not calculated					
Ni	231.604†	-0.7	-0.3594	ug/L	0.45384	-0.3594	ug/L	0.45384	126.27%
	QC value within limits for Ni	231.604	Recovery =	Not calculated					
Pb	220.353†	0.9	2.065	ug/L	1.7342	2.065	ug/L	1.7342	83.97%
	QC value within limits for Pb	220.353	Recovery =	Not calculated					
Sb	206.836†	11.2	43.27	ug/L	4.254	43.27	ug/L	4.254	9.83%
	QC value within limits for Sb	206.836	Recovery =	Not calculated					
Se	196.026†	0.9	10.77	ug/L	12.690	10.77	ug/L	12.690	117.79%
	QC value within limits for Se	196.026	Recovery =	Not calculated					
SiO2	251.603†	9.7	2.990	ug/L	1.9741	2.990	ug/L	1.9741	66.01%
	QC value within limits for SiO2	251.603	Recovery =	Not calculated					
Sr	421.552†	71.7	0.047	ug/L	0.0052	0.047	ug/L	0.0052	11.02%
	QC value within limits for Sr	421.552	Recovery =	Not calculated					
Ti	334.940†	62.5	0.296						

Zn 206.200† 0.6 0.560 ug/L 1.1158 0.560 ug/L 1.1158 199.29%
QC value within limits for Zn 206.200 Recovery = Not calculated
All analyte(s) passed QC.

Sequence No.: 34
 Sample ID: C101104-39
 Analyst: Walker
 Initial Sample Wt:
 Dilution:

Autosampler Location: 46
 Date Collected: 11/29/2010 9:23:06 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: C101104-39

Analyte Back Pressure Flow
 All 115.0 kPa 0.80 L/min

Mean Data: C101104-39

Analyte	Mean Corrected Intensity	Calib Conc. Units	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2457038.2	100.3 %	0.69			0.69%
Sc Radial	308236.7	102.8 %	0.30			0.29%
Ag 328.068†	-7.7	0.6251 ug/L	0.29626	0.6251 ug/L	0.29626	47.39%
Al 396.153†	124.1	12.78 ug/L	4.028	12.78 ug/L	4.028	31.51%
As 193.696†	-8.5	-87.93 ug/L	20.398	-87.93 ug/L	20.398	23.20%
Ba 233.527†	706.9	58.24 ug/L	0.647	58.24 ug/L	0.647	1.11%
Be 313.107†	-12.4	-0.1542 ug/L	0.02012	-0.1542 ug/L	0.02012	13.05%
B 249.677†	-52.9	-8.125 ug/L	0.8874	-8.125 ug/L	0.8874	10.92%
Ca 317.933†	178181.5	61820 ug/L	119.5	61820 ug/L	119.5	0.19%
Cd 214.440†	1.2	0.4011 ug/L	0.13881	0.4011 ug/L	0.13881	34.61%
Co 228.616†	2.0	0.6937 ug/L	1.00068	0.6937 ug/L	1.00068	144.25%
Cr 267.716†	-3.4	0.1950 ug/L	0.40862	0.1950 ug/L	0.40862	209.52%
Cu 324.752†	150.8	1.696 ug/L	0.1872	1.696 ug/L	0.1872	11.04%
Fe 238.204†	3.6	19.26 ug/L	1.174	19.26 ug/L	1.174	6.09%
K 766.490†	1348.4	758.2 ug/L	8.08	758.2 ug/L	8.08	1.07%
Mg 285.213†	48168.5	7501 ug/L	19.1	7501 ug/L	19.1	0.26%
Mn 257.610†	842.9	5.854 ug/L	0.0996	5.854 ug/L	0.0996	1.70%
Mo 202.031†	4.6	8.515 ug/L	2.0989	8.515 ug/L	2.0989	24.65%
Na 589.592†	12971.5	2374 ug/L	6.5	2374 ug/L	6.5	0.27%
Ni 231.604†	-0.7	-0.5673 ug/L	0.39511	-0.5673 ug/L	0.39511	69.64%
Pb 220.353†	-6.8	-14.45 ug/L	15.762	-14.45 ug/L	15.762	109.08%
Sb 206.836†	0.7	0.3580 ug/L	7.87577	0.3580 ug/L	7.87577	>999.9%
Se 196.026†	8.3	93.69 ug/L	89.308	93.69 ug/L	89.308	95.32%
SiO2 251.603†	21549.7	6501 ug/L	36.8	6501 ug/L	36.8	0.57%
Sr 421.552†	558363.4	365.7 ug/L	0.61	365.7 ug/L	0.61	0.17%
Ti 334.940†	-30.9	-0.146 ug/L	0.0009	-0.146 ug/L	0.0009	0.65%
Tl 190.801†	-1.4	-8.914 ug/L	18.2000	-8.914 ug/L	18.2000	204.17%
V 290.880†	-91.4	-4.855 ug/L	0.9480	-4.855 ug/L	0.9480	19.53%
Zn 206.200†	11.1	10.10 ug/L	0.905	10.10 ug/L	0.905	8.97%

Sequence No.: 35
Sample ID: C101104-42
Analyst: Walker
Initial Sample Wt:
Dilution:

Autosampler Location: 47
Date Collected: 11/29/2010 9:26:10 AM
Data Type: Original
Initial Sample Vol:
Sample Prep Vol:

Nebulizer Parameters: C101104-42

Analyte	Back Pressure	Flow
All	116.0 kPa	0.80 L/min

Mean Data: C101104-42

Analyte	Mean Corrected Intensity	Conc. Units	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2470553.0	100.8 %	0.65			0.65%
Sc Radial	305263.7	101.8 %	0.73			0.72%
Ag 328.068†	7.2	0.8197 ug/L	0.10840	0.8197 ug/L	0.10840	13.22%
Al 396.153†	118.5	12.36 ug/L	2.141	12.36 ug/L	2.141	17.33%
As 193.696†	-8.5	-87.11 ug/L	12.105	-87.11 ug/L	12.105	13.90%
Ba 233.527†	657.8	54.18 ug/L	0.518	54.18 ug/L	0.518	0.96%
Be 313.107†	-4.3	-0.1180 ug/L	0.03207	-0.1180 ug/L	0.03207	27.18%
B 249.677†	-97.6	-14.98 ug/L	1.074	-14.98 ug/L	1.074	7.17%
Ca 317.933†	166955.9	57930 ug/L	101.5	57930 ug/L	101.5	0.18%
Cd 214.440†	-0.9	-0.2387 ug/L	0.42905	-0.2387 ug/L	0.42905	179.74%
Co 228.616†	3.9	1.132 ug/L	0.7406	1.132 ug/L	0.7406	65.40%
Cr 267.716†	-0.2	0.6502 ug/L	0.46015	0.6502 ug/L	0.46015	70.77%
Cu 324.752†	150.6	1.690 ug/L	0.1011	1.690 ug/L	0.1011	5.98%
Fe 238.204†	1.7	6.165 ug/L	20.2141	6.165 ug/L	20.2141	327.90%
K 766.490†	1251.8	705.2 ug/L	12.55	705.2 ug/L	12.55	1.78%
Mg 285.213†	45627.0	7106 ug/L	7.8	7106 ug/L	7.8	0.11%
Mn 257.610†	558.3	3.793 ug/L	0.0638	3.793 ug/L	0.0638	1.68%
Mo 202.031†	3.5	6.343 ug/L	3.1646	6.343 ug/L	3.1646	49.89%
Na 589.592†	12300.1	2251 ug/L	7.1	2251 ug/L	7.1	0.32%
Ni 231.604†	0.5	0.1280 ug/L	1.53321	0.1280 ug/L	1.53321	>999.9%
Pb 220.353†	-5.8	-12.31 ug/L	6.183	-12.31 ug/L	6.183	50.21%
Sb 206.836†	-0.4	-3.977 ug/L	15.9780	-3.977 ug/L	15.9780	401.78%
Se 196.026†	8.1	91.94 ug/L	7.294	91.94 ug/L	7.294	7.93%
SiO2 251.603†	21295.7	6425 ug/L	51.9	6425 ug/L	51.9	0.81%
Sr 421.552†	521992.1	341.9 ug/L	0.56	341.9 ug/L	0.56	0.16%
Ti 334.940†	-36.5	-0.172 ug/L	0.0104	-0.172 ug/L	0.0104	6.05%
Tl 190.801†	-2.8	-14.69 ug/L	8.742	-14.69 ug/L	8.742	59.49%
V 290.880†	-95.7	-4.914 ug/L	0.6231	-4.914 ug/L	0.6231	12.68%
Zn 206.200†	9.1	8.301 ug/L	1.5808	8.301 ug/L	1.5808	19.04%

Sequence No.: 36
 Sample ID: C101104-45
 Analyst: Walker
 Initial Sample Wt:
 Dilution:

Autosampler Location: 48
 Date Collected: 11/29/2010 9:29:15 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: C101104-45

Analyte	Back Pressure	Flow
All	117.0 kPa	0.80 L/min

Mean Data: C101104-45

Analyte	Mean Corrected Intensity	Conc. Units	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2460841.8	100.4 %	0.32			0.31%
Sc Radial	306743.6	102.3 %	0.47			0.46%
Ag 328.068†	-5.6	0.8195 ug/L	0.10558	0.8195 ug/L	0.10558	12.88%
Al 396.153†	144.7	17.71 ug/L	1.189	17.71 ug/L	1.189	6.71%
As 193.696†	-11.0	-113.1 ug/L	18.48	-113.1 ug/L	18.48	16.34%
Ba 233.527†	680.9	55.94 ug/L	0.444	55.94 ug/L	0.444	0.79%
Be 313.107†	-10.7	-0.1608 ug/L	0.04452	-0.1608 ug/L	0.04452	27.68%
B 249.677†	-139.5	-21.42 ug/L	0.748	-21.42 ug/L	0.748	3.49%
Ca 317.933†	193745.3	67220 ug/L	148.0	67220 ug/L	148.0	0.22%
Cd 214.440†	-0.5	-0.1099 ug/L	0.58757	-0.1099 ug/L	0.58757	534.47%
Co 228.616†	4.8	1.376 ug/L	0.3281	1.376 ug/L	0.3281	23.84%
Cr 267.716†	-2.7	0.3770 ug/L	0.40958	0.3770 ug/L	0.40958	108.64%
Cu 324.752†	179.5	2.037 ug/L	0.0680	2.037 ug/L	0.0680	3.34%
Fe 238.204†	7.0	42.60 ug/L	11.253	42.60 ug/L	11.253	26.41%
K 766.490†	1406.6	783.2 ug/L	22.97	783.2 ug/L	22.97	2.93%
Mg 285.213†	51228.5	7978 ug/L	34.9	7978 ug/L	34.9	0.44%
Mn 257.610†	7573.8	54.96 ug/L	0.313	54.96 ug/L	0.313	0.57%
Mo 202.031†	4.7	8.836 ug/L	2.4791	8.836 ug/L	2.4791	28.06%
Na 589.592†	14319.4	2620 ug/L	21.6	2620 ug/L	21.6	0.82%
Ni 231.604†	0.8	0.3026 ug/L	1.34288	0.3026 ug/L	1.34288	443.84%
Pb 220.353†	-9.5	-20.72 ug/L	5.638	-20.72 ug/L	5.638	27.21%
Sb 206.836†	-0.4	-4.274 ug/L	9.2926	-4.274 ug/L	9.2926	217.41%
Se 196.026†	7.9	89.32 ug/L	17.986	89.32 ug/L	17.986	20.14%
SiO2 251.603†	22036.8	6646 ug/L	38.0	6646 ug/L	38.0	0.57%
Sr 421.552†	726174.7	475.7 ug/L	0.80	475.7 ug/L	0.80	0.17%
Ti 334.940†	-43.1	-0.204 ug/L	0.0380	-0.204 ug/L	0.0380	18.64%
Tl 190.801†	-5.4	-25.83 ug/L	9.371	-25.83 ug/L	9.371	36.29%
V 290.880†	-83.2	-4.851 ug/L	0.0098	-4.851 ug/L	0.0098	0.20%
Zn 206.200†	10.7	9.636 ug/L	0.3237	9.636 ug/L	0.3237	3.36%

Sequence No.: 37
 Sample ID: C101104-48
 Analyst: Walker
 Initial Sample Wt:
 Dilution:

Autosampler Location: 49
 Date Collected: 11/29/2010 9:32:21 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: C101104-48

Analyte Back Pressure Flow
 All 115.0 kPa 0.80 L/min

Mean Data: C101104-48

Analyte	Mean Corrected Intensity	Calib Conc. Units	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2406777.6	98.20 %	0.378			0.39%
Sc Radial	304138.7	101.5 %	0.44			0.43%
Ag 328.068†	-29.8	0.2610 ug/L	0.12669	0.2610 ug/L	0.12669	48.54%
Al 396.153†	432.7	22.72 ug/L	7.673	22.72 ug/L	7.673	33.78%
As 193.696†	-6.7	-71.83 ug/L	35.786	-71.83 ug/L	35.786	49.82%
Ba 233.527†	182.8	13.32 ug/L	0.340	13.32 ug/L	0.340	2.55%
Be 313.107†	-21.9	-0.4090 ug/L	0.07679	-0.4090 ug/L	0.07679	18.78%
B 249.677†	-291.9	-44.82 ug/L	0.421	-44.82 ug/L	0.421	0.94%
Ca 317.933†	664781.3	230700 ug/L	624.8	230700 ug/L	624.8	0.27%
Cd 214.440†	13.9	4.376 ug/L	1.5634	4.376 ug/L	1.5634	35.73%
Co 228.616†	3.8	1.559 ug/L	1.7806	1.559 ug/L	1.7806	114.22%
Cr 267.716†	-7.8	0.7133 ug/L	0.16727	0.7133 ug/L	0.16727	23.45%
Cu 324.752†	191.9	1.574 ug/L	0.1746	1.574 ug/L	0.1746	11.09%
Fe 238.204†	3.4	9.994 ug/L	21.4882	9.994 ug/L	21.4882	215.01%
K 766.490†	2742.0	1588 ug/L	19.6	1588 ug/L	19.6	1.23%
Mg 285.213†	138330.4	21540 ug/L	45.7	21540 ug/L	45.7	0.21%
Mn 257.610†	358.2	1.935 ug/L	0.0423	1.935 ug/L	0.0423	2.19%
Mo 202.031†	5.1	6.364 ug/L	4.5777	6.364 ug/L	4.5777	71.93%
Na 589.592†	57405.0	10530 ug/L	35.0	10530 ug/L	35.0	0.33%
Ni 231.604†	0.5	-1.361 ug/L	1.3409	-1.361 ug/L	1.3409	98.53%
Pb 220.353†	-13.0	-21.10 ug/L	5.368	-21.10 ug/L	5.368	25.44%
Sb 206.836†	-0.3	-2.375 ug/L	5.5627	-2.375 ug/L	5.5627	234.24%
Se 196.026†	17.4	187.6 ug/L	30.13	187.6 ug/L	30.13	16.06%
SiO2 251.603†	35665.4	10790 ug/L	45.1	10790 ug/L	45.1	0.42%
Sr 421.552†	Saturated2					
Ti 334.940†	-60.6	-0.287 ug/L	0.1106	-0.287 ug/L	0.1106	38.58%
Tl 190.801†	-1.5	-12.83 ug/L	12.590	-12.83 ug/L	12.590	98.13%
V 290.880†	-24.3	-4.517 ug/L	0.4754	-4.517 ug/L	0.4754	10.53%
Zn 206.200†	969.5	918.2 ug/L	6.14	918.2 ug/L	6.14	0.67%

Sequence No.: 38
 Sample ID: C101104-51
 Analyst: Walker
 Initial Sample Wt:
 Dilution:

Autosampler Location: 50
 Date Collected: 11/29/2010 9:36:27 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: C101104-51

Analyte Back Pressure Flow
 All 115.0 kPa 0.80 L/min

Mean Data: C101104-51

Analyte	Mean Corrected Intensity	Calib Conc. Units	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2408282.9	98.26 %	0.611			0.62%
Sc Radial	307518.9	102.6 %	0.31			0.30%
Ag 328.068†	-5.0	0.4163 ug/L	0.27957	0.4163 ug/L	0.27957	67.15%
Al 396.153†	863.0	109.9 ug/L	7.51	109.9 ug/L	7.51	6.84%
As 193.696†	-4.0	-42.69 ug/L	44.792	-42.69 ug/L	44.792	104.92%
Ba 233.527†	236.6	17.58 ug/L	0.129	17.58 ug/L	0.129	0.73%
Be 313.107†	130.6	0.1228 ug/L	0.06077	0.1228 ug/L	0.06077	49.50%
B 249.677†	-243.3	-37.36 ug/L	0.625	-37.36 ug/L	0.625	1.67%
Ca 317.933†	695940.9	241500 ug/L	1440.7	241500 ug/L	1440.7	0.60%
Cd 214.440†	56.1	17.27 ug/L	1.063	17.27 ug/L	1.063	6.16%
Co 228.616†	14.7	4.067 ug/L	0.7363	4.067 ug/L	0.7363	18.10%
Cr 267.716†	-8.2	-0.4839 ug/L	0.25201	-0.4839 ug/L	0.25201	52.08%
Cu 324.752†	1794.1	18.41 ug/L	0.167	18.41 ug/L	0.167	0.91%
Fe 238.204†	399.2	2834 ug/L	26.3	2834 ug/L	26.3	0.93%
K 766.490†	3004.7	1751 ug/L	27.1	1751 ug/L	27.1	1.55%
Mg 285.213†	130982.0	20400 ug/L	91.0	20400 ug/L	91.0	0.45%
Mn 257.610†	289637.1	2114 ug/L	2.7	2114 ug/L	2.7	0.13%
Mo 202.031†	4.0	4.232 ug/L	4.9058	4.232 ug/L	4.9058	115.91%
Na 589.592†	60609.7	11120 ug/L	54.8	11120 ug/L	54.8	0.49%
Ni 231.604†	12.3	5.084 ug/L	1.3503	5.084 ug/L	1.3503	26.56%
Pb 220.353†	-12.7	-20.80 ug/L	9.285	-20.80 ug/L	9.285	44.63%
Sb 206.836†	-3.9	-17.56 ug/L	12.268	-17.56 ug/L	12.268	69.87%
Se 196.026†	12.6	128.5 ug/L	45.04	128.5 ug/L	45.04	35.06%
SiO2 251.603†	55261.2	16700 ug/L	134.9	16700 ug/L	134.9	0.81%
Sr 421.552†	Saturated2					
Ti 334.940†	-67.6	-0.320 ug/L	0.0610	-0.320 ug/L	0.0610	19.08%
Tl 190.801†	-1.7	-18.58 ug/L	5.070	-18.58 ug/L	5.070	27.29%
V 290.880†	-10.8	-4.991 ug/L	0.2314	-4.991 ug/L	0.2314	4.64%
Zn 206.200†	3782.2	3576 ug/L	10.1	3576 ug/L	10.1	0.28%

Sequence No.: 39
 Sample ID: 1011110-BLK1
 Analyst: Walker
 Initial Sample Wt:
 Dilution:

Autosampler Location: 51
 Date Collected: 11/29/2010 9:40:32 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: 1011110-BLK1

Analyte	Back Pressure	Flow
All	115.0 kPa	0.80 L/min

Mean Data: 1011110-BLK1

Analyte	Mean Corrected Intensity	Conc. Units	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2461510.6	100.4 %	0.70			0.70%
Sc Radial	305779.0	102.0 %	0.26			0.26%
Ag 328.068†	11.2	0.1876 ug/L	0.16094	0.1876 ug/L	0.16094	85.77%
Al 396.153†	4.5	0.8627 ug/L	0.32655	0.8627 ug/L	0.32655	37.85%
As 193.696†	-1.7	-17.28 ug/L	16.154	-17.28 ug/L	16.154	93.50%
Ba 233.527†	-0.1	-0.0158 ug/L	0.06803	-0.0158 ug/L	0.06803	431.37%
Be 313.107†	-13.7	-0.0502 ug/L	0.02168	-0.0502 ug/L	0.02168	43.24%
B 249.677†	1.7	0.2604 ug/L	0.20551	0.2604 ug/L	0.20551	78.93%
Ca 317.933†	18.2	6.292 ug/L	1.2048	6.292 ug/L	1.2048	19.15%
Cd 214.440†	1.4	0.4344 ug/L	0.36928	0.4344 ug/L	0.36928	85.01%
Co 228.616†	0.3	0.0801 ug/L	0.35490	0.0801 ug/L	0.35490	443.18%
Cr 267.716†	-1.4	-0.2238 ug/L	0.18988	-0.2238 ug/L	0.18988	84.84%
Cu 324.752†	131.5	1.373 ug/L	0.2056	1.373 ug/L	0.2056	14.97%
Fe 238.204†	3.6	25.34 ug/L	27.226	25.34 ug/L	27.226	107.44%
K 766.490†	53.2	30.36 ug/L	19.775	30.36 ug/L	19.775	65.13%
Mg 285.213†	-2.4	-0.3762 ug/L	0.25427	-0.3762 ug/L	0.25427	67.59%
Mn 257.610†	-0.4	-0.0055 ug/L	0.00392	-0.0055 ug/L	0.00392	71.12%
Mo 202.031†	1.5	3.013 ug/L	4.5586	3.013 ug/L	4.5586	151.28%
Na 589.592†	16.1	2.880 ug/L	5.3295	2.880 ug/L	5.3295	185.05%
Ni 231.604†	0.4	0.2480 ug/L	0.55624	0.2480 ug/L	0.55624	224.29%
Pb 220.353†	-0.9	-2.038 ug/L	5.2686	-2.038 ug/L	5.2686	258.57%
Sb 206.836†	0.9	3.521 ug/L	10.6168	3.521 ug/L	10.6168	301.57%
Se 196.026†	2.1	25.20 ug/L	20.186	25.20 ug/L	20.186	80.12%
SiO2 251.603†	-20.0	-5.989 ug/L	1.3464	-5.989 ug/L	1.3464	22.48%
Sr 421.552†	190.6	0.125 ug/L	0.0329	0.125 ug/L	0.0329	26.34%
Ti 334.940†	6.6	0.031 ug/L	0.0101	0.031 ug/L	0.0101	32.38%
Tl 190.801†	-1.7	-7.168 ug/L	10.5085	-7.168 ug/L	10.5085	146.61%
V 290.880†	-92.6	-3.264 ug/L	0.3441	-3.264 ug/L	0.3441	10.54%
Zn 206.200†	1.0	0.913 ug/L	0.7842	0.913 ug/L	0.7842	85.93%

Sequence No.: 40
 Sample ID: 1011110-BS1
 Analyst: Walker
 Initial Sample Wt:
 Dilution:

Autosampler Location: 52
 Date Collected: 11/29/2010 9:43:37 AM
 Data Type: Original
 Initial Sample Vol:
 Sample Prep Vol:

Nebulizer Parameters: 1011110-BS1

Analyte	Back Pressure	Flow
All	115.0 kPa	0.80 L/min

Mean Data: 1011110-BS1

Analyte	Mean Corrected Intensity	Conc. Units	Std.Dev.	Sample Conc. Units	Std.Dev.	RSD
Sc Axial	2487294.4	101.5 %	1.16			1.15%
Sc Radial	308446.0	102.9 %	1.36			1.32%
Ag 328.068†	6097.0	101.0 ug/L	1.93	101.0 ug/L	1.93	1.91%
Al 396.153†	49187.6	10570 ug/L	57.6	10570 ug/L	57.6	0.55%
As 193.696†	7.3	85.64 ug/L	23.649	85.64 ug/L	23.649	27.61%
Ba 233.527†	1202.2	99.47 ug/L	1.563	99.47 ug/L	1.563	1.57%
Be 313.107†	28276.2	103.1 ug/L	1.30	103.1 ug/L	1.30	1.26%
B 249.677†	-75.1	-11.54 ug/L	0.767	-11.54 ug/L	0.767	6.65%
Ca 317.933†	29864.0	10340 ug/L	147.3	10340 ug/L	147.3	1.42%
Cd 214.440†	338.8	103.4 ug/L	1.66	103.4 ug/L	1.66	1.60%
Co 228.616†	428.1	100.7 ug/L	1.88	100.7 ug/L	1.88	1.86%
Cr 267.716†	609.8	97.34 ug/L	1.211	97.34 ug/L	1.211	1.24%
Cu 324.752†	9448.2	98.77 ug/L	1.682	98.77 ug/L	1.682	1.70%
Fe 238.204†	1458.0	10400 ug/L	140.0	10400 ug/L	140.0	1.35%
K 766.490†	18476.7	10440 ug/L	52.8	10440 ug/L	52.8	0.51%
Mg 285.213†	67456.4	10510 ug/L	34.7	10510 ug/L	34.7	0.33%
Mn 257.610†	13769.3	100.1 ug/L	1.00	100.1 ug/L	1.00	1.00%
Mo 202.031†	48.1	98.20 ug/L	0.896	98.20 ug/L	0.896	0.91%
Na 589.592†	57746.1	10580 ug/L	48.4	10580 ug/L	48.4	0.46%
Ni 231.604†	177.5	98.22 ug/L	1.543	98.22 ug/L	1.543	1.57%
Pb 220.353†	42.5	95.05 ug/L	5.502	95.05 ug/L	5.502	5.79%
Sb 206.836†	19.1	70.64 ug/L	12.987	70.64 ug/L	12.987	18.39%
Se 196.026†	44.4	525.5 ug/L	22.69	525.5 ug/L	22.69	4.32%
SiO2 251.603†	34.6	-8.079 ug/L	2.1092	-8.079 ug/L	2.1092	26.11%
Sr 421.552†	839968.0	550.2 ug/L	0.97	550.2 ug/L	0.97	0.18%
Ti 334.940†	-60.2	-0.285 ug/L	0.0276	-0.285 ug/L	0.0276	9.68%
Tl 190.801†	29.4	117.7 ug/L	3.36	117.7 ug/L	3.36	2.86%
V 290.880†	2639.0	90.64 ug/L	2.664	90.64 ug/L	2.664	2.94%
Zn 206.200†	109.3	101.3 ug/L	2.46	101.3 ug/L	2.46	2.43%

Matrix Recovery Check: 1011110-BS1

Analyte	Expected Conc.	Measured Conc.	Std. Dev.	Units	Recovery (%)
Al 396.153	10100	10570	57.590	ug/L	104.6
Ca 317.933	10110	10340	147.280	ug/L	102.4
Fe 238.204	10130	10400	139.962	ug/L	102.7
K 766.490	10130	10440	52.807	ug/L	103.1
Mg 285.213	10100	10510	34.657	ug/L	104.0
Na 589.592	10100	10580	48.422	ug/L	104.7
Ag 328.068	100.2	101.0	1.933	ug/L	100.8
As 193.696	82.72	85.64	23.649	ug/L	102.9
Ba 233.527	99.98	99.47	1.563	ug/L	99.5
Be 313.107	99.95	103.1	1.299	ug/L	103.2
Cd 214.440	100.4	103.4	1.657	ug/L	103.0
Co 228.616	100.1	100.7	1.877	ug/L	100.6
Cr 267.716	99.78	97.34	1.211	ug/L	97.6
Cu 324.752	101.4	98.77	1.682	ug/L	97.4
Mn 257.610	99.99	100.1	1.001	ug/L	100.1
Mo 202.031	103.0	98.20	0.896	ug/L	95.2
Ni 231.604	100.2	98.22	1.543	ug/L	98.0
Pb 220.353	97.96	95.05	5.502	ug/L	97.1
Sb 206.836	103.5	70.64	12.987	ug/L	67.1
Se 196.026	525.2	525.5	22.690	ug/L	100.1