



CB&I
2700 Chandler Avenue, Building C
Las Vegas, NV 89120
Tel: +1 702 795 0515
Fax: +1 702 795 8210
www.CBI.com

May 9, 2014

Dania Zinner
USEPA; Region 8
1595 Wynkoop Street (8EPR-SR)
Denver, CO 80202-1129

Document ID #: 3019-05122014-1

Dear Ms. Zinner:

EPA CONTRACT NUMBER EP-W-10-033
TASK ORDER NUMBER 3019
QA SUPPORT FOR THE LIBBY ASBESTOS SITE

Enclosed please find the Summary Asbestos On-site Audit Report for the on-site audit performed on March 25-26, 2014 at EMSL Analytical, Inc. in Cinnaminson, New Jersey. This report and the accompanying checklist are deliverables under Task 5 of Task Order 3019.

If you have any questions, please feel free to contact me.

Sincerely,

A handwritten signature in black ink, appearing to read 'Timothy L. Vonnahme', written in a cursive style.

Timothy L. Vonnahme
Audit Group Manager, QATS Program
CB&I Federal Services, LLC
Phone: (702) 895-8729
E-Mail Address: timothy.vonnahme@cbifederalservices.com

cc: Administrative Contracting Officer (letter only)
Audit Group Files



The Quality Assurance Technical Support (QATS) contract is operated by Shaw Environmental, Inc., a CB&I Company.
The QATS Program's Quality Management System is certified to the ISO 9001:2008 International Standard

REPORT
FOR
TASK ORDER NUMBER 3019
QUALITY ASSURANCE SUPPORT FOR THE LIBBY ASBESTOS SITE
SUMMARY ASBESTOS ON-SITE AUDIT REPORT

EMSL Analytical, Inc. (Cinnaminson, NJ)

Prepared by:

**The Data Auditing Group
Quality Assurance Technical Support Program
CB&I Federal Services LLC
2700 Chandler Avenue
Las Vegas, Nevada 89120**

April 24, 2014

QATS Contract Number: EP-W-10-033

Prepared for:

**Dania Zinner
Task Order Manager**

**Region 8
U.S. Environmental Protection Agency
1595 Wynkoop Street
Denver, CO 80202**

TABLE OF CONTENTS

LABORATORY INFORMATION AND AUDIT SCOPE 3

EXECUTIVE SUMMARY 4

AUDIT FINDINGS..... 5

 Sample Receipt, Storage, Log-in, and Chain-of-Custody (COC)..... 5

 Phase Contrast Microscopy (PCM) 5

 Indirect and Direct Preparation of Air Filter and Dust Samples for TEM Analysis 5

 Transmission Electron Microscopy (TEM) Analysis..... 6

 Polarized Light Microscopy (PLM) Analysis..... 6

 Data Management 7

 Quality Control and Quality Assurance 7

CORRECTIVE ACTION APPLIED FROM THE PREVIOUS AUDIT FINDINGS 7

CONCLUSIONS 9

ATTACHMENT

 Libby-Specific Asbestos Laboratory On-site Audit Checklist (EPA only)

LABORATORY INFORMATION AND AUDIT SCOPE

This report summarizes the findings of an Asbestos on-site laboratory audit of the EMSL Analytical, Inc. Laboratory in Cinnaminson, New Jersey conducted on March 25-26, 2014. The audit was conducted in support of the United States Environmental Protection Agency (EPA) to assess the performance of laboratories supporting Libby Superfund Site activities. CB&I Federal Services LLC Quality Assurance Technical Support (QATS) staff participation in the on-site audit and subsequent preparation of this report was performed under Task 5, Task Order 3019, QATS Contract EP-W-10-033.

Detailed information regarding the subject laboratory is as follows:

Date of On-site:	March 25-26, 2014
Laboratory:	EMSL Analytical, Inc. 200 Route 130 North Cinnaminson, NJ 08077 800.220.3675
Special Projects Manager:	Robyn Denton
Audit Team	
US EPA:	Charlie Appleby, ASB CLP Project Manager
CB&I QATS:	Michael P. Lenkauskas, CQA, Lead Auditor

The Audit Team, which was comprised of EPA Analytical Services Branch (ASB) and CB&I Federal Services LLC QATS personnel, performed the technical and evidentiary aspects of the on-site audit. The technical part of the audit involved an evaluation of the laboratory's facilities, personnel, and capabilities to process samples and data as described in the Libby-specific guidance documents. Processes evaluated included sample receipt, sample storage, sample tracking, sample preparation, sample analysis, data review, and data package assembly. Laboratory instrumentation and equipment were inspected for proper maintenance and calibration, and laboratory personnel were interviewed to determine proficiency in their assigned responsibilities. Specific instrumentation and areas inspected included sample receiving, Phase Contrast Microscopy (PCM), Polarized Light Microscopy (PLM), Transmission Electron Microscopy (TEM), and the laboratory's capability to provide the required hardcopy and electronic data deliverables (EDDs).

The evidentiary part of the evaluation involved an assessment of laboratory documentation for accuracy, completeness, and defensibility. The Laboratory Quality Assurance Manual (QAM) and Standard Operating Procedures (SOPs) were assessed for availability and accuracy to observed procedures, and instrument calibration and maintenance logbooks were reviewed for completeness, traceability, and accuracy. During the course of the audit, the Libby-Specific Asbestos Laboratory On-site Audit Checklist was completed by the QATS Auditor. The checklist is provided as an attachment to this report (EPA only).

EXECUTIVE SUMMARY

An Asbestos on-site audit of EMSL Analytical, Inc. in Cinnaminson, New Jersey was performed on March 25-26, 2014 in support of Libby Superfund Site activities. Areas assessed included facilities, equipment, personnel, and documentation as related to the laboratory's capability to process samples collected from the Libby Superfund Site and to test for Asbestos and other fibers in accordance with Libby-specific requirements. The on-site audit identified six (6) deficiencies which are summarized below by laboratory area:

Sample Receipt, Storage, Log-in, and Chain-of-Custody – Written procedures describing remote login procedures are not available.

Phase Contrast Microscopy – The results of daily reference slide analyses are not documented; therefore, results failing acceptance criteria are not being recorded.

Indirect and Direct Preparation of Air Filter and Dust Samples – Analyst training for preparing various media (i.e., bark, duff, and FBAS) for TEM analyses has not been performed. (Note: this is a repeat defect). Balances do not contain updated calibration labels as required in the project and laboratory procedures.

Data Management – The electronic spreadsheet used to track electronic and hardcopy deliverable due dates did not have a column for tracking the due dates of hardcopy deliverables. (Note: this was corrected prior to the close of the audit.)

Quality Control and Quality Assurance – The Internal Audit SOP was not reviewed and revised at the 3-year frequency required in the Laboratory's QAM.

With the exception of the deficiencies noted above and in the following report, the on-site evaluation revealed that the EMSL Analytical, Inc. laboratory in Cinnaminson, NJ to have sufficient facilities, equipment, and staff to effectively analyze samples in accordance with the Libby-specified methodologies. All staff and management were cooperative, readily answered questions by the Audit Team, and appeared to be responsive to the identified audit findings.

AUDIT FINDINGS

Sample Receipt, Storage, Log-in, and Chain-of-Custody (COC)

The sample receipt area, located in the reception area, was clean and well organized. The Audit Team interviewed the Special Projects Data Coordinator on the procedures used to inspect, process, and login samples collected from the Libby Superfund Site. The Special Projects Data Coordinator demonstrated a clear understanding of the process for sample inspection, processing, and distribution. One deficiency related to the lack of an SOP for performing remote login of samples received at the EMSL facilities in Denver, CO and New York, NY was identified:

1. Written procedures describing how the Special Projects Data Coordinator remotely logs in and generates the necessary documents for samples collected at the Libby Superfund Site and distributed to the EMSL Analytical Laboratories in Denver, CO and New York, NY are not available. Although Section 5.4.2 of the Laboratory's Sample Chain-of-Custody references a Remote Login SOP, the SOP was not available. The requirement that instructions or procedures for the activities affecting the quality of analytical services be developed by management is described in Section 5.4.1 of the Laboratory's QAM. (Checklist No. 4.6.1)

Recommended Corrective Action – Ensure that written procedures for the remote login procedures used to process samples and generate documents for the EMSL Analytical Laboratories in Denver, CO and New York, NY are available.

Phase Contrast Microscopy (PCM)

The PCM area was clean and organized, with adequate equipment and instrumentation for preparing and analyzing air samples by PCM. The analyst demonstrated a clear understanding of the applicable techniques for inspecting and preparing air filters samples as described in the applicable Libby-specific guidance documents. One deficiency concerning quality control analyses was identified:

2. The results from a daily reference slide, which is required to be read by each analyst prior to analyzing samples, are entered into an electronic spreadsheet which indicates whether or not the result is within criteria and, therefore, whether the analyst can continue to analyze client samples. However, there is no mechanism to document results that fail criteria. As a result, there is no way to determine whether the analyst has entered multiple numbers until one is entered that meets criteria or analyzed client samples after multiple failed attempts. The requirements that the analysis fall within the acceptable limits before analysis may proceed, and that a Corrective Action Response (CAR) is initiated for those results that fall outside the limits are described in Section 8.5.3 of the Laboratory's Asbestos and Other Fibers by PCM SOP. (Checklist No. 5.8.1)

Recommended Corrective Action – Ensure that results of failed daily reference analyses are recorded in a permanent manner and that a CAR is initiated.

Indirect and Direct Preparation of Air Filter and Dust Samples for TEM Analysis

The TEM preparation area was clean and organized. Adequate equipment and instrumentation were available for preparing air, dust, water, tree bark, and duff samples for TEM analysis using the appropriate direct and indirect preparation techniques. The analyst demonstrated a clear understanding of the applicable techniques for inspecting and preparing air filters samples as

described in the applicable Libby-specific guidance documents. Two deficiencies concerning equipment calibration and training were identified:

3. During the interview, it was determined that the analyst has received training for the direct and indirect preparation of TEM air filter samples received from the Libby Superfund Site; however, training in the preparation of other media that could be received during the upcoming sampling season (i.e., bark, duff, and FBAS) had not been received. The training requirements for laboratory personnel are described in Section 4.2.3.3 of the Site-wide Quality Assurance Reference Document (QARD, Rev. 1), and Section 5.2.2 of the Laboratory's Quality Assurance Manual (QAM). (Checklist No. 10.2.1)

Recommended Corrective Action – Ensure that additional personnel are trained in the preparation of all media that could be received from the Libby Superfund site for TEM analysis.

4. At the time of the audit, the labels on the analytical balances indicated that the balances were past due the 12 month recalibration by an outside vendor. However, it was later determined that the labels applied were incorrect. Supporting documentation indicated that the balances had been calibrated as required within the last 12 months. The requirement that all balances be labeled with the date of the certification, initials of the individual performing the calibration and certification, and the date the next service is to be performed are described in Section 12.2.1 of the project-specific SOP SRC-Libby-01 (Rev. 3) and Section 5.5.3 of the Laboratory's QAM. (Checklist Nos. 6.4.4, 6.15.1, 8.4.4.5, and 8.16.1)

Note: This finding also applies to the balances used to weigh samples for the PLM-GRAV and PLM-VE procedures.

Recommended Corrective Action – Ensure that all balances have calibration labels with a sticker indicating the correct date of the certification, initials of the individual performing the calibration and certification, and the date the next service is to be performed.

Transmission Electron Microscopy (TEM) Analysis

The area was clean and well organized. The TEM instruments used to support the project were well-maintained, calibrated at the specified frequencies, and equipped with digital photography capabilities. The TEM analyst interviewed demonstrated a clear understanding of the applicable techniques for identifying and recording structures as described in the applicable Libby-specific guidance documents. No deficiencies concerning TEM analyses were identified.

Polarized Light Microscopy (PLM) Analysis

The PLM area has three work stations, each equipped to analyze samples received from the Libby Superfund Site. Each work station is equipped with a stereomicroscope, functional HEPA hood, polarized light microscope, refractive index (RI) liquids, and tools for manipulating samples. The PLM area was clean and organized; the instrumentation well-maintained; and the quality of the documentation acceptable. The analyst interviewed demonstrated a clear understanding of PLM instrument maintenance and calibration, and sample preparation, analysis, and documentation.

As a follow-up to the recent PLM inter-laboratory study, the Audit Team asked that the analyst reanalyze samples which were reported as weakly discordant from the original analysis. From the reanalysis, a confirmed result was reported. No PLM deficiencies were identified.

Data Management

Data management activities associated with the analysis of samples collected at the Libby Superfund Site are performed by EMSL's Special Projects group, which reviews all records of sample receipt, preparation, and analysis for accuracy, compliance, and completeness. This group is also responsible for generating the hardcopy and electronic deliverables for special projects, including the Libby Superfund Site, for all participating EMSL Analytical laboratories. The Special Projects area was clean and well organized, and the procedures to ensure data completeness and integrity adequate. The Special Projects Data Coordinator responsible for data management activities clearly described her duties with respect to data review and the generation of data deliverables. One deficiency concerning the tracking of hardcopy deliverable due dates was identified:

5. The electronic spreadsheet used by the Special Projects Data Coordinator to track both electronic and hardcopy deliverable due dates did not have a column for tracking the due dates of hardcopy deliverables. Having this column would allow the laboratory the ability to notify the client when data would not be submitted by the specified due date. The requirement to provide the customer with information regarding deliverables that cannot be provided on time is described in Section 4.7.4 of the Laboratory's QAM. (Checklist No. 9.2.4.2)

Recommended Corrective Action – Prior to the audit debriefing, a column indicating hardcopy deliverable due dates was added to the deliverables tracking spreadsheet. Therefore, no further corrective action is necessary.

Quality Control and Quality Assurance

The Audit Team interviewed the Quality Assurance Officer (QAO), reviewed the laboratory's QAM and SOPs, and performed a cursory review of the laboratory's air monitoring results, non-conformance reports, laboratory certifications, internal audit reports, and the training files of select laboratory personnel. The QAO demonstrated an understanding of and commitment to the laboratory's current quality system. One deficiency concerning the timely review of laboratory SOPs was identified:

6. The Internal Audit SOP has not been revised since 2010, which exceeds the laboratory SOP review cycle requirement of a minimum of every three years. The requirement to review controlled documents once every three years to determine their continued suitability is described in Section 4.3.1.5 of the laboratory's QAM. (Checklist No. 10.3.1.1)

Recommended Corrective Action – With the exception of the QAM, which is reviewed on an annual basis, ensure that all controlled documents are reviewed at a minimum of every three years to determine the continued suitability.

CORRECTIVE ACTION APPLIED FROM THE PREVIOUS AUDIT FINDINGS

The on-site laboratory evaluation included an assessment of the findings reported in the previous Summary Asbestos On-site Audit Report for the on-site audit performed on October 8, 2013. Of the three findings identified in the previous on-site audit, the laboratory has completely

addressed two (66.7%) and partially addressed one (33.3%). The following are the findings identified during the previous on-site audit, the laboratory's verbatim responses to the findings (where applicable), and observations made during the current on-site audit.

Indirect and Direct Preparation of Air Filter and Dust Samples

1. With the departure of a key staff member, the laboratory no longer has adequate personnel with the training necessary to prepare duff and tree bark samples for analysis by TEM. The training requirements for laboratory personnel are described in Section 4.2.3.3 of the Site-wide Quality Assurance Reference Document (QARD, Rev. 0) and Section 5.2.2 of the laboratory's Quality Assurance Manual (QAM).

Recommended Corrective Action – Ensure that properly trained personnel are available to prepare tree bark, duff, and other samples received from the Libby Superfund site.

EMSL Corrective Action Response (12/20/2013): *The lab has begun to train additional staff in Libby sample preparation. Currently, three staff members are trained in Libby specific preparation of TEM filters. These staff members are: Robyn Denton, Leslie McCluskey- Eissing and Kim Ford. As different media is received, media specific training will be completed.*

Effectiveness Check (03/25-26/2014): This deficiency has been partially addressed. Although some training has been completed for the indirect preparation of air and dust samples, training for other media received from the Libby Superfund Site (i.e., bark, duff, and FBAS indirect preparation) has not yet been performed.

Polarized Light Microscopy (PLM) Analysis

2. One of the PLM microscopes used to analyze Libby samples incorrectly utilized a 530 nm compensator plate, rather than a 550 nm compensator plate as specified in the Libby project-specific procedure. Although this deviation is recorded on the bench sheet, it is not described in sufficient detail and is not described in the applicable data package narratives. The requirement to use a 550 nm compensator plate is described in Section 10.3.1.12 of both the PLM-VE (SRC-Libby-03, Rev. 3) and PLM-Grav (SRCLibby-01, Rev. 3).

Recommended Corrective Action – Ensure all deviations from project-specific requirements are described in sufficient detail in the applicable data package narratives.

EMSL Corrective Action Response (12/20/2013): *The lab has ensured that any deviations from project specific requirements are described in the sample case narrative.*

Effectiveness Check (03/25-26/2014): This deficiency has been completely addressed. Since this finding was identified, the project-specific SOPs have been modified through Laboratory Modification LB-000097 to allow for 530-550 nm compensator plates, which is consistent with the NVLAP requirements.

Data Management

3. The laboratory was not including the data package completeness checklists with the data deliverables. These checklists are provided with both the PLM-VE and PLM-GRAV EDD templates. The requirement to provide a completed checklist with each PLM-VE,

PLM-GRAV, and NIOSH 9002 hardcopy (scanned) data deliverable is described in the "Data Pkg Checklist" tab of each of the applicable EDD templates.

Recommended Corrective Action – Ensure that data package checklists are provided with each PLM-VE, PLM-GRAV, and NIOSH 9002 hardcopy deliverable.

EMSL Corrective Action Response (12/20/2013): *Robyn Denton met with staff members regarding the use of the data package checklists. Since the previous on-site audit, EMSL has been including the checklist with all PLM-VE, PLM_GRAV and NIOSH 9002 hard copy deliverables. Please see attachment 3A, which is from a recent PLM VE Data package.*

Effectiveness Check (03/25-26/2014): This deficiency has been completely addressed.

CONCLUSIONS

An Asbestos on-site audit was performed at EMSL Analytical, Inc. in Cinnaminson, New Jersey on March 25-26, 2014. The audit involved an assessment of the laboratory's facility, instrumentation, personnel, and laboratory procedures to process samples received from the Libby Superfund Site in Libby, Montana. The on-site audit identified the following six (6) deficiencies:

- Written procedures describing remote login procedures are not available.
- The results of daily PCM reference slide analyses are not documented.
- Additional training is needed for preparing various media for TEM analyses. (Note: this is a partial repeat defect).
- Balances do not contain updated calibration labels as required in the project and laboratory procedures.
- The electronic spreadsheet used to track electronic and hardcopy deliverable due dates did not have a column for tracking the due dates of hardcopy deliverables. (Note: this was corrected prior to the audit close).
- The Internal Audit SOP was not reviewed at the 3-year frequency required in the Laboratory's QAM.

With the exception of the deficiencies noted above, the on-site evaluation revealed the laboratory to have sufficient facilities, equipment, and staff to effectively analyze samples in accordance with the specified methodologies and Libby-specific protocol. All staff and management were cooperative, readily answered questions by the Audit Team, and appeared to be responsive to the identified audit findings.

ATTACHMENT

Libby-Specific Asbestos Laboratory On-site Audit Checklist (EPA Only)

LIBBY-SPECIFIC ASBESTOS LABORATORY ON-SITE AUDIT CHECKLISTMethod: Not ApplicableDate(s) of On-site: March 25-26, 2014Laboratory: EMSL Analytical, Inc.Address: 200 Route 130 North
Cinnaminson, NJ 08077Telephone: (800) 220-3675Laboratory Personnel Contacted

Name	Title
Robyn Denton	Special Projects Manager/PCM Analyst
Charles LaCerra	Special Projects/Sample Receiving Manager
Garret Vliet	PLM Supervisor
Meghan Smollock	Special Projects Data Coordinator
Melissa Klinedinst	PLM QC Group Leader
Leslie McCluskey-Eissing	TEM Analyst

Evaluation Team

Name	Title
Charlie Appleby	EPA-ASB, CLP Project Manager
Michael P. Lenkauskas	CB&I Federal Services LLC, Lead Auditor

LIBBY-SPECIFIC ASBESTOS LABORATORY ON-SITE AUDIT CHECKLIST

Method: Not Applicable

Date(s) of On-site: March 25-26, 2014

TABLE OF CONTENTS

1.0 LABORATORY STATUS & CAPABILITIES..... 1

2.0 LABORATORY SECURITY..... 1

3.0 PROJECT INITIATION/PROJECT MANAGEMENT 1

4.0 SAMPLE RECEIPT, LOG-IN, STORAGE, & TRACKING 2

5.0 PHASE CONTRAST MICROSCOPY (PCM) 4

6.0 TRANSMISSION ELECTRON MICROSCOPY (TEM) GRID PREPARATION..... 6

7.0 TEM ANALYSIS 14

8.0 POLARIZED LIGHT MICROSCOPY (PLM)..... 19

9.0 DATA MANAGEMENT 28

10.0 QUALITY ASSURANCE/QUALITY CONTROL 30

LIBBY-SPECIFIC ASBESTOS LABORATORY ON-SITE AUDIT CHECKLIST

Method: Not ApplicableDate(s) of On-site: March 25-26, 2014

1.0 LABORATORY STATUS & CAPABILITIES		Yes	No	Comments
1.1 Which of the following capabilities does the laboratory possess:				
1.1.1	Phase Contrast Microscopy (PCM)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	Full Service Lab.
1.1.2	Polarized Light Microscopy (PLM)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
1.1.3	Transmission Electron Microscopy (TEM)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
1.1.4	Others (list)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
1.2 Is the laboratory currently receiving samples from Libby Superfund Site Operable Units?		<input checked="" type="checkbox"/>	<input type="checkbox"/>	
If "YES," complete the following table:				
Method	Media	Comments		
ISO 10312	Various	Air, tree bark, duff, water, FBAS, etc...		
AHERA	Air			
ASTM	Dust			
PLM	Soil			
PCM	Air			
2.0 LABORATORY SECURITY		Yes	No	Comments
2.1 Are visitors required to sign in?		<input checked="" type="checkbox"/>	<input type="checkbox"/>	
2.2 Are all entrances to the laboratory secure?		<input checked="" type="checkbox"/>	<input type="checkbox"/>	
3.0 PROJECT INITIATION/PROJECT MANAGEMENT		Yes	No	Comments
3.1 Are there designated project managers or a project management team to ensure received samples are properly processed?		<input checked="" type="checkbox"/>	<input type="checkbox"/>	Robyn Denton
3.2 Are project-specific requirements and procedures communicated to laboratory personnel and available for reference:				All project specific documents are available in the CDM Smith eRoom, and laboratory documents are available on the laboratory's server.
3.2.1	Project-specific SOPs?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
3.2.2	Laboratory Modifications?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
3.2.3	SAP Analytical Summaries?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
3.2.4	Project-specific Electronic Data Deliverables (EDDs)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
3.2.5	Other (list)?	<input type="checkbox"/>	<input checked="" type="checkbox"/>	
Additional Comments:				

LIBBY-SPECIFIC ASBESTOS LABORATORY ON-SITE AUDIT CHECKLIST

Method: Not ApplicableDate(s) of On-site: March 25-26, 2014

4.0 SAMPLE RECEIPT, LOG-IN, STORAGE, & TRACKING		Yes	No	Comments
4.1	Is the sample receiving area adequate, clean, and orderly?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Personnel Interviewed				
Name		Title		Experience
Meghan Smollock		Special Projects Data Coordinator		2 Years
Charles LaCerra		Sample Receiving Manager		14 Years
Robyn Denton		Special Projects Manager		13 years
4.2 Sample Receipt		Yes	No	Comments
4.2.1	Is there a sample custodian and designated alternate responsible for sample receipt and log-in?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	See Additional Comments below.
4.2.2	Is the sample custodian or alternate available to receive and log-in samples at any time delivery services are operating?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
4.2.3	Are sample shipping containers opened in a HEPA hood (as necessary) to both minimize personal exposure and safeguard against laboratory contamination?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
4.2.4	Does the sample custodian verify and record the following when inspecting shipments and reviewing documentation:			
4.2.4.1	Presence and condition of custody seals?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
4.2.4.2	The SAP analytical summary is referenced or provided?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
4.2.4.3	Presence or absence of Chain-of-Custody (COC) records?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
4.2.4.4	Presence or absence of air bill sticker(s)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
4.2.4.5	Sample condition?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
4.2.4.6	Presence of packaging or packing material which could compromise samples (i.e., vermiculite & polystyrene)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
4.2.4.7	Problems/discrepancies between samples, documentation, client requests, etc.?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
4.2.4.8	Bulk and air samples received separately?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
4.2.5	Are COC records signed and dated at the time of sample receipt?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	Time stamped upon arrival.
4.2.6	Is a system in place to contact the client in case of absent documentation or discrepancies between COCs, client requests, etc.?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
4.2.7	Are subsequent resolutions to problems and discrepancies documented?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	As described in project-specific COC SOP.
4.3 Sample Identification				
4.3.1	Are sample receipt identification logbooks, or a LIMS, used to log-in samples and assign unique laboratory identification numbers?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	LIMS is used.
4.3.1.1	Does the logbook or logging system serve as a direct cross-reference between laboratory ID numbers and client ID numbers?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Additional Comments:				
Meghan is responsible for inspection, review, login and distribution of all samples received from the Libby Superfund Site. In addition, Meghan also remotely logs in Libby samples received at the EMSL locations in Denver, CO and New York, NY.				

LIBBY-SPECIFIC ASBESTOS LABORATORY ON-SITE AUDIT CHECKLIST

Method: Not ApplicableDate(s) of On-site: March 25-26, 2014

4.0 SAMPLE RECEIPT, LOG-IN, STORAGE, & TRACKING	Yes	No	Comments
4.4 Sample Storage			
4.4.1 Are storage facilities sufficient?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	See Additional Comments below.
4.4.2 Is the sample storage area secured to prevent entry of unauthorized personnel?	<input type="checkbox"/>	<input type="checkbox"/>	
4.4.3 Is a logbook or other means used to record sample locations?	<input type="checkbox"/>	<input type="checkbox"/>	
4.4.4 Are samples easy to locate from logbook references? Select and find a previously analyzed sample (Sample no. <u>N/A</u>)	<input checked="" type="checkbox"/>	<input type="checkbox"/>	Did not track an individual sample but inspected off-site facility as described below.
4.5 Sample Tracking			
4.5.1 Is a system in place to keep track of samples entering and leaving the storage, sample preparation, and analysis areas?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
4.5.2 Is the retention and/or disposal of unused portions of samples and prepared samples documented?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
4.5.2.1 Are project-specific retention and/or disposal requirements communicated and followed?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
4.6 Standard Operating Procedures (SOPs)			
4.6.1 Are the applicable laboratory SOPs available and followed by laboratory personnel (list)?	<input type="checkbox"/>	<input checked="" type="checkbox"/>	Refer to Finding No. 1 of the On-site Audit Report.
Document Title	Control No.	Description	
Sample Chain-of-Custody	Rev. 3 (8/16/2012)		
4.7 Document Control:	Yes	No	Comments
4.7.1 Are all logbooks, notebooks, forms, or other laboratory documents legible, accurate, and complete (list)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Document Title	Description/Comments		
Libby-specific QC Log books	Used to assign QC samples for Libby projects		
Additional Comments:			
<p>Samples are archived for long-term and short-term storage at a remote facility, which the Audit Team visited on the afternoon of the first day of the audit. The storage facility, which is also used to archive data, is secure and organized with EMSL personnel on-site during normal business hours.</p>			

LIBBY-SPECIFIC ASBESTOS LABORATORY ON-SITE AUDIT CHECKLIST

Method: Not ApplicableDate(s) of On-site: March 25-26, 2014

5.0 PHASE CONTRAST MICROSCOPY (PCM)		Yes	No	Comments
5.1	Does the laboratory perform PCM analyses on samples received from the Libby Superfund site? If "NO," proceed to Section 6.0 of the checklist.	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
5.2	Is the PCM area adequate, clean, and orderly?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
5.3	Are steps taken to prevent the cross-contamination of equipment, supplies, and reagents?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Personnel Interviewed				
Name		Title		Experience
Robyn Denton		Special Projects Manager		13 years
5.4 Methods and Guidance Documents		Yes	No	Comments
5.4.1	Are the applicable guidance documents available for reference:			
5.4.1.1	NIOSH Method 7400 (Issue 2), 1994?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
5.4.1.2	Other (list)?	<input type="checkbox"/>	<input checked="" type="checkbox"/>	
5.4.2	Are project-specific requirements communicated to laboratory personnel and available for reference:			
5.4.2.1	Most recent revision of Laboratory Modification LB-000015?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	All project specific documents and EDD templates are available in the CDM Smith eRoom.
5.4.2.2	SOP EPA-Libby-08?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
5.4.2.3	SAP Analytical Summaries?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
5.4.2.4	Project-specific Electronic Data Deliverables (EDDs)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
5.4.2.5	Other (list)?	<input type="checkbox"/>	<input checked="" type="checkbox"/>	
5.5 Equipment				
5.5.1	Ventilation Hoods:			
5.5.1.1	Checked routinely and recorded in a permanent logbook?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
5.5.2	Are the microscopes used to analyze samples equipped with the following:			
5.5.2.1	Positive phase contrast, with green or blue filter?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
5.5.2.2	Adjustable field iris?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
5.5.2.3	Eyepiece (8 to 10X)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
5.5.2.4	Phase magnification (40 to 45X)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
5.5.2.5	Walton-Beckett Graticule?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
5.5.2.6	Stage micrometer with 0.01 mm subdivisions?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
5.5.3	Are microscope and phase ring alignment checks conducted daily?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
5.5.4	Is resolution periodically checked using an HSE/NPL slide?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
5.5.5	Are maintenance and calibration activities recorded in microscope-specific logbooks?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Additional Comments:				

LIBBY-SPECIFIC ASBESTOS LABORATORY ON-SITE AUDIT CHECKLIST

Method: Not ApplicableDate(s) of On-site: March 25-26, 2014

5.0 PHASE CONTRAST MICROSCOPY (PCM)	Yes	No	Comments
5.6 Sample Preparation			
5.6.1 Are filters prepared as described in the applicable method(s)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
5.6.2 Are filters that are visibly overloaded (>25%) or that contain loose debris prepared indirectly as described in SOP EPA-Libby-08?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
5.7 Sample Analysis			
5.7.1 Are the appropriate counting rules used (A or B)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	"A" rules are followed.
5.7.2 How are the fields and fibers tracked and recorded? <u>Calibrated double counters</u>			
5.8 Quality Control			
5.8.1 Is each analyst provided a minimum of one reference slide per work day?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	Refer to Finding No. 2 of the On-site Audit Report.
5.8.2 Are recounts analyzed at a frequency of 1 per 10 samples analyzed?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
5.8.2.1 For count pairs not within acceptance limits, are associated samples recounted?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
5.9 Standard Operating Procedures (SOPs)			
5.9.1 Are the applicable laboratory SOPs available and followed by laboratory personnel (list)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Document Title	Control No.	Description	
Asbestos & Other Fibers by PCM	Rev. 14.7 (6/14/2013)		
5.10 Document Control	Yes	No	Comments
5.10.1 Are all logbooks, notebooks, forms, or other laboratory documents legible, accurate, and complete (list)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Document Title	Description/Comments		
PCM Calibration Logbook	Track microscope and counter calibrations		
Additional Comments:			

LIBBY-SPECIFIC ASBESTOS LABORATORY ON-SITE AUDIT CHECKLIST

Method: Not ApplicableDate(s) of On-site: March 25-26, 2014

6.0 TRANSMISSION ELECTRON MICROSCOPY (TEM) GRID PREPARATION	Yes	No	Comments
6.1 Are the grid preparation areas adequate, clean, and orderly?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.2 Are bulk samples prepared in an area separate from that used to prepare air and dust samples?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.3 Are steps taken to prevent the cross-contamination of equipment, supplies, and reagents?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Personnel Interviewed			
Name	Title		Experience
Leslie McCluskey-Eissing	TEM Analyst		2 years
Robyn Denton	Special Projects Manager		13 Years
6.4 Equipment & Supplies	Yes	No	Comments
6.4.1 Ventilation Hoods:			
6.4.1.1 Checked routinely and recorded in a permanent logbook?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.4.2 Drying oven:			
6.4.2.1 Checked routinely and recorded in a permanent logbook?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	Two separate ovens calibrated to 60 and 80 degrees Celsius.
6.4.3 Muffle furnace:			
6.4.3.1 Checked routinely and recorded in a permanent logbook?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	Four ovens calibrated to 480 degrees Celsius.
6.4.4 Analytical balances:			
6.4.4.1 Checked routinely and recorded in a permanent logbook?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	Refer to Finding No. 4 of the On-site Audit Report.
6.4.4.2 Calibrated within the last 12 months by a certified technician?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.4.5 Plasma Asher (refer to the most recent revision of Laboratory Modification LB-000085):			
6.4.5.1 Calibrated at least quarterly and recorded in a permanent logbook?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.4.6 Sputter Coater (Vacuum evaporator):			
6.4.6.1 Checked routinely and recorded in a permanent logbook?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.4.7 Filtration Apparatus (for indirect preparation):			
6.4.7.1 Are disposable funnels used (record catalogue #)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	The EFA of each bag is checked and recorded on the bag with the most recent EFA identified as having an EFA of 364.9 mm ² .
6.4.7.2 Has the Effective Filtration Area (EFA) been determined and recorded for each funnel lot?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.4.8 TEM Grids:			
6.4.8.1 Is documentation for average grid opening determination available?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Additional Comments:			

LIBBY-SPECIFIC ASBESTOS LABORATORY ON-SITE AUDIT CHECKLIST

Method: Not ApplicableDate(s) of On-site: March 25-26, 2014

6.0 TRANSMISSION ELECTRON MICROSCOPY (TEM) GRID PREPARATION	Yes	No	Comments
6.5 Direct and Indirect Preparation Methodology			
6.5.1 What method(s) does the laboratory use to prepare air and dust samples for TEM analysis:			
6.5.1.1 40 CFR, Chapter 1, Part 763, Subpart E - AHERA?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.5.1.2 ISO 10312:1195 E - Determination of Asbestos Fibers?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.5.1.3 ASTM D 5755-09 - Micro vacuum Sampling and Indirect Analysis of Dust by TEM?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.5.1.4 Others (list)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	EPA 100.2
6.5.2 Are project-specific requirements communicated to laboratory personnel and available for reference:			
6.5.2.1 Laboratory Modifications?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.5.2.2 Project-specific SOPs?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.5.2.3 SAP Analytical Summaries?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.5.2.4 Other (list)?	<input type="checkbox"/>	<input checked="" type="checkbox"/>	All project specific documents are available in the CDM Smith eRoom.
6.6 Sample Inspection			
6.6.1 Are air filter cassettes carefully wet-wiped prior to being transferred to the clean preparation area for inspection?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.6.2 Are air filter samples which are visibly overloaded, exhibit uneven loading, or contain loose debris, prepared indirectly?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.6.3 Are all ambient air samples dried upon receipt at the laboratory prior to preparation and analysis (refer to the most recent revision of Laboratory Modification LB-000055B)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.6.3.1 Is a drying blank (DB) prepared?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.7 Direct Preparation of MCE and Polycarbonate (PC) Filters			
6.7.1 Are MCE filters collapsed using either a Di-Methyl Formamide (DMF) or Acetone Atmosphere (AA) technique (describe technique)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Refer to Laboratory Modification LB-000091 for acetone use.			
6.7.2 Is plasma etching performed on collapsed MCE filters?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.7.2.1 Is a 5 to 10% layer of the collapsed surface removed during etching?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	5% is etched.
6.7.3 Are collapsed MCE filters and secured polycarbonate filters transferred to a vacuum evaporator for carbon coating?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.7.4 Are excised filter sections placed on the appropriately labeled TEM grids and cleared using a Jaffe Washer or an equivalent technique (describe)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.7.5 Are samples checked for remaining filter residue after clearing?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.7.5.1 If residue remains, is condensation washing or an equivalent technique used (describe technique)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Additional Comments:			

LIBBY-SPECIFIC ASBESTOS LABORATORY ON-SITE AUDIT CHECKLIST

Method: Not ApplicableDate(s) of On-site: March 25-26, 2014

6.0 TRANSMISSION ELECTRON MICROSCOPY (TEM) GRID PREPARATION	Yes	No	Comments
6.8 Indirect Sample Preparation of Air and Dust Samples			
6.8.1 Are the applicable Libby guidance documents available for reference:			
6.8.1.1 SOP EPA-Libby-08 - Indirect Preparation of Air and Dust Sample for TEM Analysis?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.8.1.2 Most recent revision of Laboratory Modification LB-000091?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.8.2 Sample filtration:			
6.8.2.1 Are the applicable SAP Analytical Summaries reviewed to determine whether or not filter samples must be ashed?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.8.2.2 Are cassettes examined for loose material?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.8.2.2.1 If loose material or uneven loading is not evident, is a portion of the air samples retained?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.8.2.2.2 If loose material is evident, is the loose material filtered along with the air filter?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.8.2.3 Ashing (if applicable):			
6.8.2.3.1 Are filters covered with aluminum foil and placed in a plasma asher?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.8.2.3.2 Is the plasma asher operated at minimum power?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.8.2.3.3 Is 100% ashing confirmed by visual observation?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.8.2.4 Are air filters, loose material, dust, or ash, rinsed into a beaker and brought to a final volume of 100 mL with particle-free water?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.8.2.4.1 Adjusted to a pH of 3-4 with a 10% solution of Glacial Acetic Acid (refer to Laboratory Modification LB-000091 exemption)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.8.2.4.2 Sonicated for 3 minutes and allowed to settle for 2 minutes prior to filtering?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.8.2.5 Are the appropriate aliquots of filtrate passed through a <u>disposable</u> 25 mm filter assembly with a 0.2 µm MCE filter with a 5.0 µm MCE support pad?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.8.2.6 Is a secondary filter loading of between 10% and 25% achieved?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.8.3 Are serial dilutions performed as necessary?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.8.4 Are all dilution volumes recorded on an indirect preparation bench sheet and provided in the associated data deliverable?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.8.5 Are TEM grids prepared as described in Section 6.7 of this checklist?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Additional Comments:			

LIBBY-SPECIFIC ASBESTOS LABORATORY ON-SITE AUDIT CHECKLIST

Method: Not ApplicableDate(s) of On-site: March 25-26, 2014

6.0 TRANSMISSION ELECTRON MICROSCOPY (TEM) GRID PREPARATION	Yes	No	Comments
6.9 Water Sample Preparation			
6.9.1 What method(s) does the laboratory use to prepare water samples for TEM analysis:			
6.9.1.1 EPA Method 100.2 - Determination of Asbestos Structures Over 10 µm in Length in Drinking Water?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.9.1.2 EPA Method 100.1 - Determination of Asbestos Fibers Drinking Water?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.9.1.3 Others (describe)? _____	<input type="checkbox"/>	<input checked="" type="checkbox"/>	
6.9.2 Is sample preparation performed in accordance with the most recent revision of Laboratory Modification LB-000020:			UV light source is broken. However, a new one is on order and will arrive before sampling begins.
6.9.2.1 Do samples undergo treatment with ozone/UV light?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.9.2.2 Are Sample aliquots of no less than 1 mL poured through MCE or PC filters with a pore size of 0.22 µm or smaller?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.9.2.3 Are the following processes QC'd by a second person and documented by an initial and date on the preparation bench sheet:			
6.9.2.3.1 Ozone treatment?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.9.2.3.2 Filtration?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.9.2.3.3 Assignment to Petri dishes?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.9.2.3.4 Placement on glass slides for etching and carbon coating?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.9.2.3.5 Grid preparation?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.9.2.4 After aliquots have been filtered, is the remaining volume archived in its original container until the Laboratory Controller (LC) requests it be filtered for permanent archival?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.9.3 Are all dilution volumes recorded on an indirect preparation bench sheet, and provided in the associated data deliverable?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.9.4 Are TEM grids prepared as described in Section 6.7 of this checklist?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Additional Comments:			

LIBBY-SPECIFIC ASBESTOS LABORATORY ON-SITE AUDIT CHECKLIST

Method: Not ApplicableDate(s) of On-site: March 25-26, 2014

6.0 TRANSMISSION ELECTRON MICROSCOPY (TEM) GRID PREPARATION	Yes	No	Comments
6.10 OU3 Tree Bark Sample Preparation			
6.10.1 Are the applicable Libby guidance documents available for reference:			
6.10.1.1 EPA-Libby-2012-12 - Sampling and Analysis of Tree Bark for Asbestos?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.10.2 Drying and Ashing:			
6.10.2.1 Are the diameter and thickness of the tree bark samples measured and recorded to an accuracy of ± 2 mm?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.10.2.2 Is the entire tree bark sample, which can sometimes be multiple core samples, weighed and placed in an oven for drying?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.10.2.2.1 Is the sample dried at 80°C until the weight stabilizes (a minimum of 6 hours) and weighed?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.10.2.3 Is the tree bark sample then covered, placed in a muffle furnace at 450°C for 18 hours (or until all organic matter has been removed), and weighed?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.10.2.3.1 Is the furnace ramped from 0°C to 450°C?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.10.3 Acid Treatment:			
6.10.3.1 After adding approximately 1-2 mL of DI water, is 10-20 mL of concentrated HCL added until no further reaction is visible (approx. 3-5 minutes)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.10.3.2 Are samples diluted, transferred to a 100 mL container (with lid), and brought to a final volume of 100 mL with fiber-free DI water?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.10.3.3 Are samples capped, inverted 5-6 times, and sonicated for 2 minutes in preparation for filtering?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.10.4 Filtration:			
6.10.4.1 Are 5-20 mL of solution transferred to a second container and brought to a volume of 100 mL with fiber-free DI water?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.10.4.2 Are dilutions agitated (inverted 5-6 times) and filtered through a 47 mm MCE filter (0.45 μ m pore size)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.10.4.2.1 Are additional dilutions prepared if the loading on the filter appears either too heavy (>20%) or too light?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.10.5 Are all dilution volumes recorded on an indirect preparation bench sheet, and provided in the associated data deliverable?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.10.6 Are TEM grids prepared as described in Section 6.7 of this checklist?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Additional Comments:			

LIBBY-SPECIFIC ASBESTOS LABORATORY ON-SITE AUDIT CHECKLIST

Method: Not ApplicableDate(s) of On-site: March 25-26, 2014

6.0 TRANSMISSION ELECTRON MICROSCOPY (TEM) GRID PREPARATION	Yes	No	Comments
6.11 OU3 Duff Sample Preparation			
6.11.1 Are the applicable Libby guidance documents available for reference:			
6.11.1.1 EPA-Libby-2012-11 - Sampling and Analysis of Duff for Asbestos?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.11.2 Drying and Ashing:			
6.11.2.1 Are the appropriate number of aluminum trays weighed and tared?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.11.2.1.1 For tracking purposes, is each tray marked with a unique number?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.11.2.2 Are trays filled to approximately $\frac{3}{4}$, dried at 60°C until the weight stabilizes a minimum of 10 hours, and weighed?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.11.2.3 Are dried duff samples transferred to covered pans and placed in a muffle furnace at 450°C for 18 hours, or until all organic matter has been removed, and weighed?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.11.2.4 Are ashed samples transferred to Zip-lock bags and homogenized?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.11.2.4.1 If an individual sample was split between multiple trays, was it combined into one Zip-lock bag?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.11.3 Acid Treatment:			
6.11.3.1 After adding approximately 1-2 mL of DI water to 0.25 grams (measured to ± 0.01 g) of ashed sample, is 10-20 mL of concentrated HCL added until no further reaction is visible (approx. 3-5 minutes)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.11.3.2 Are samples diluted, transferred to a 100 mL container (with lid) and brought to a final volume of 100 mL with fiber-free DI water?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.11.3.3 Are samples capped, inverted 5-6 times, and sonicated for 2 minutes in preparation for filtering?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.11.4 Filtration:			
6.11.4.1 Is 0.1-1.0 mL of solution transferred to a second container and brought to a volume of 100 mL with fiber-free DI water?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.11.4.2 Are dilutions agitated (inverted 5-6 times) and filtered through a 47 mm MCE filter (0.45 μ m pore size)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.11.4.2.1 Are additional dilutions prepared if the loading on the filter appears either too heavy (>20%) or too light?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.11.5 Are all dilution volumes recorded on an indirect preparation bench sheet, and provided in the associated data deliverable?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.11.6 Are TEM grids prepared as described in Section 6.7 of this checklist?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Additional Comments:			

LIBBY-SPECIFIC ASBESTOS LABORATORY ON-SITE AUDIT CHECKLIST

Method: Not ApplicableDate(s) of On-site: March 25-26, 2014

6.0 TRANSMISSION ELECTRON MICROSCOPY (TEM) GRID PREPARATION		Yes	No	Comments
6.14 Quality Control Samples				
6.14.1 Are quality control samples prepared at the described frequency:				
6.14.1.1	Are laboratory blanks (LBs) prepared at a frequency of 4% or with each preparation batch, whichever is more frequent?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.14.1.2	Are re-preparations prepared at a frequency of 1%?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.14.2 Is there a system in place to track and assign quality control analyses for samples associated with SAP Summaries that have different frequency requirements than those found in the most recent revision of Laboratory Modification LB-000029?		<input checked="" type="checkbox"/>	<input type="checkbox"/>	
6.15 Standard Operating Procedures (SOPs)				
6.15.1 Are the applicable laboratory SOPs available and followed by laboratory personnel (list)?		<input type="checkbox"/>	<input checked="" type="checkbox"/>	Refer to Finding No. 4 of the On-site Audit Report.
Document Title		Control No.		Description
Carbon Coater		Rev. 4.1 (5/2/2012)		
Direct Transfer Filter Prep		Rev. 1.4 (8/8/2013)		
TEM GO Measurement		Rev. 5 (3/22/2013)		
Plasma Asher		Rev. 4.1 (12/11/2013)		
6.16 Document Control		Yes	No	Comments
6.16.1 Are all logbooks, notebooks, forms, or other laboratory documents legible, accurate, and complete (list)?		<input type="checkbox"/>	<input checked="" type="checkbox"/>	
Document Title		Description/Comments		
Balance calibration		Balances are calibrated using three sets of NIST traceable weights.		
Muffle Furnace Calibration		Calibrated at the specified frequencies.		
Drying Furnace Calibration		Calibrated at the specified frequencies.		
Additional Comments:				

LIBBY-SPECIFIC ASBESTOS LABORATORY ON-SITE AUDIT CHECKLIST

Method: Not ApplicableDate(s) of On-site: March 25-26, 2014

7.0 TEM ANALYSIS		Yes	No	Comments
7.1	Are TEM areas adequate, clean, and orderly?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.2	Are steps taken to prevent the cross-contamination of equipment, supplies, and reagents?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Personnel Interviewed				
Name		Title		Experience
Leslie McCluskey-Eissing		TEM Analyst		4 years
Robyn Denton		Special Projects Manager		13 Years
7.3 Methods and Guidance Documents		Yes	No	Comments
7.3.1	What method(s) does the laboratory use to analyze samples by TEM:			
7.3.1.1	40 CFR, Chapter 1, Part 763, Subpart E (AHERA)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.3.1.2	ISO 10312:1995 E - Determination of Asbestos Fibers?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.3.1.3	ASTM D 5755-09 - Microvacuum Sampling and Indirect Analysis of Dust by TEM?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.3.1.4	EPA Method 100.2 - Determination of Asbestos Structures Over 10 µm in Length in Drinking Water?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.3.1.5	Others (list)? _____	<input type="checkbox"/>	<input checked="" type="checkbox"/>	
7.3.2	Are project-specific requirements communicated to laboratory personnel and available for reference:			
7.3.2.1	Laboratory Modifications?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.3.2.2	Project-specific SOPs?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.3.2.3	SAP Analytical Summaries?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.3.2.4	Project-specific Electronic Data Deliverables (EDDs)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.3.2.5	Other (list)? _____	<input type="checkbox"/>	<input checked="" type="checkbox"/>	All project specific documents and EDD templates are available in the CDM Smith eRoom.
7.4 TEM Instrumentation				
7.4.1	Does TEM instrumentation meet the following requirements:			
7.4.1.1	Capable of being operated at between 80 and 120 kV?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.4.1.2	Electron diffraction (ED) and energy dispersive X-ray (EDX) capabilities?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.4.1.3	Fluorescent screen with an inscribed or overlaid calibrated scale?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.4.2	Are the instruments equipped with thin film or Beryllium windows (list below)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.4.3	Are all routine and non-routine maintenance activities recorded in instrument-specific logbooks?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Instrument No.	Make	Model	Capabilities	
04-1	100CX-2	JOEL	Horizontal detector w/Be window	
04-3	1200EX	JOEL	Horizontal detector w/Be window	
04-5	100CX-2	JOEL	Horizontal detector w/Be window	
04-6	1200EX2	JOEL	Horizontal detector w/Be window	
Additional Comments:				

LIBBY-SPECIFIC ASBESTOS LABORATORY ON-SITE AUDIT CHECKLIST

Method: Not ApplicableDate(s) of On-site: March 25-26, 2014

7.0 TEM ANALYSIS	Yes	No	Comments
7.5 Instrument Calibration (Laboratory Modification LB-00085)			
7.5.1 Is microscope alignment performed <u>daily</u> :			
7.5.1.1 Centering of electron beam?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.5.1.2 Electron beam is properly stigmated on either side of crossover?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.5.1.3 Image properly focused?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.5.2 Is the TEM screen magnification calibrated <u>monthly</u> ?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.5.3 Is the camera constant calibrated <u>monthly</u> ?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.5.4 Is the spot size diameter determined to be less than 250 nm <u>quarterly</u> ?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.5.5 Is the low beam dose (≥ 15 seconds for Chrysotile) verified <u>quarterly</u> ?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.5.6 EDXA System:			
7.5.6.1 Is X-ray energy versus channel for two peaks (i.e., Cu/Al) checked <u>daily</u> ?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.5.6.2 Is detector resolution (Mn) checked <u>quarterly</u> ?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.5.6.3 Are K-factors performed <u>quarterly</u> :			
7.5.6.3.1 BIR-1G (Na, Mg, Al, Ca, and Fe relative to Si)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.5.6.3.2 Orthoclase (K and Al relative to Si)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.5.7 Are instrument calibration records maintained in instrument-specific logbooks?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.5.8 Are calibrations uploaded to the eRoom on a quarterly basis?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.6 Reference Materials			
7.6.1 Does the laboratory maintain a library of reference materials on Asbestos and other fiber types?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.6.2 Are instrument-specific "LA" spectra available for reference?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.7 Grid Acceptance/Rejection Criteria			
7.7.1 Grid preparation rejection criteria:			
7.7.1.1 The replica is too dark due to poor dissolution?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.7.1.2 Replica is doubled or folded?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.7.1.3 Replica has >25% obscuration rejected?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.7.1.4 Replica has <50 intact grid openings?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.7.2 Are samples associated with grids determined to be overloaded (>25%) re-prepped using the indirect-transfer technique described in SOP EPA-Libby-08 and Laboratory Modification LB-000091?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.7.2.1 For samples prepared indirectly, is a loading of between 10% and 25% achieved for the secondary filter?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Additional Comments:			

LIBBY-SPECIFIC ASBESTOS LABORATORY ON-SITE AUDIT CHECKLIST

Method: Not ApplicableDate(s) of On-site: March 25-26, 2014

7.0 TEM ANALYSIS	Yes	No	Comments
7.8 Modifications to AHERA & ASTM D5755:			
7.8.1 Most recent revision of Laboratory Modification LB-000031:			
7.8.1.1 Are structures classified as fibers (F), bundles (B), clusters (C), or matrices (M)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.8.1.2 Unless identified as a "close call," are NAMs not recorded?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.8.1.3 Is the designation "ND" used to document when no structures are detected in a grid opening?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.8.1.4 Are fibers, bundles, clusters and matrices only recorded when they contain individual constituent fibers meeting the aspect ratio criterion?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.8.1.5 The overall aspect ratio of bundles, clusters, and matrices, may have any value?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.8.1.6 Are non-countable structures recorded, but not counted, for informational purposes?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.8.1.7 Is the entire length recorded for structures originating in one grid opening and extending to an adjacent grid opening?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.8.1.8 Are the actual lengths and widths of fibers, bundles, clusters and matrices recorded?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.8.1.9 For disperse matrices and clusters, is the length of the longest protruding structure recorded?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.8.1.10 For analyses with less than 50 grid openings (GOs) counted, is selection random and are adjacent GOs avoided?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.9 Modifications to ISO Method 10312:			
7.9.1 Most recent revision of Laboratory Modification LB-000016:			
7.9.1.1 Recording of "close call" NAMS as described in the most recent revision of Laboratory Modification LB-000066?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.9.1.2 Recording of bundles only if they contain individual constituent fibers meeting the aspect ratio criteria?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.9.1.3 Recording of bundles, compact clusters, and compact matrices regardless of aspect ratio?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.9.1.4 Recording of structures that intersect countable and non-countable grid bars:			
7.9.1.4.1 Cross Grid Bar Length Doubled (XGBLD)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.9.1.4.2 Crosses Non-Countable Grid Bar Length Doubled (XNCGBLD)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.9.1.5 Recording of component structures, within non-countable structures, which do not intersect non-countable grid bars?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.9.1.6 Recording of disperse clusters and matrices?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.9.1.7 Are the recorded rules for partially obscured structures properly applied (i.e., MFO and MBO)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.9.1.8 Are the counting and recording rules for the identification of PCMe structures at "low magnification" applied?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Additional Comments:			

LIBBY-SPECIFIC ASBESTOS LABORATORY ON-SITE AUDIT CHECKLIST

Method: Not ApplicableDate(s) of On-site: March 25-26, 2014

7.0 TEM ANALYSIS	Yes	No	Comments
7.10 Common TEM Modifications:			
7.10.1 Most recent revision of Laboratory Modification LB-000066:			
7.10.1.1 Is the presence or absence of Sodium and Potassium recorded for all LA, OA, and NAM particles (NaK, NaX, XK, or XX)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.10.1.2 Is probable mineral identification code recorded for all particles?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.10.1.2.1 Are LA particles identified as WRTA, AC, TR, or AT?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.10.1.2.2 Are OA particles identified as AM, AN, CR, or NR?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.10.1.2.3 Are NAMs indicated as PY, OT, or UN?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.10.1.3 Is one SAED pattern recorded for each amphibole Asbestos type encountered per samples?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.10.1.4 Are EDS spectrum (a maximum of 5) collected for up to 5 LA, 5 NR, and 5 Close-call NAM per sample?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.10.1.5 To the extent possible, are all EDS spectra collected for sufficient time that key peaks (i.e., Na, K, and Al) can be distinguished from background?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.10.1.6 To the extent possible, are all EDS spectra collected for sufficient time that the Si peak contains 1,000 or more counts?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.10.2 Most recent revision of Laboratory Modification LB-000067:			
7.10.2.1 Is the designation "ND" used to document when no structures are detected in a grid opening?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.10.2.2 Do hardcopy bench sheets include sketches of all Asbestos structures observed, up to a maximum of 50?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.10.2.2.1 Do these sketches contain sufficient detail?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.10.2.3 Are the structure identification codes described in Tables D.1 and D.2 of ISO Method 10312 used?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.10.2.4 Are laboratory blanks assigned the sample number LQ-00001 and assigned the appropriate tag (i.e., AL1, REP1, and FBA1)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.10.2.5 Is the preparation date for blanks the date on which it is introduced to the sample train?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.10.2.6 Is the preparation date for field samples the date on which the preparation is initiated?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Additional Comments:			

LIBBY-SPECIFIC ASBESTOS LABORATORY ON-SITE AUDIT CHECKLIST

Method: Not ApplicableDate(s) of On-site: March 25-26, 2014

7.0 TEM ANALYSIS		Yes	No	Comments
7.11 Counting/Stopping rules:				
7.11.1 Are the Analytical Summaries reviewed to determine the following:				
7.11.1.1	Analytical Sensitivity?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.11.1.2	Recording rules?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.11.1.3	Stopping rules?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.11.1.4	Applicable Laboratory Modifications?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.11.1.5	Investigative or non-investigative?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.12 Quality Control Analyses (Laboratory Modification LB-000029)				
7.12.1 Are quality control samples analyzed at the required frequencies:				
7.12.1.1	Laboratory blanks - Frequency of 4%?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.12.1.2	Recount Same (RS) - Frequency of 1%?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.12.1.3	Recount Different (RD) - Frequency of 2.5%?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.12.1.4	Verified Analysis (VA) - Frequency of 1%?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.12.1.5	Re-preparations - Frequency of 1%?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.12.2 Is there a system in place to track and assign quality control analyses for samples associated with SAP Summaries that have different frequency requirements than those found in the most recent revision of Laboratory Modification LB-000029?		<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.12.3 Are samples selected for RS, RD and VA analyses in accordance with the most recent revision of Laboratory Modification LB-000029?		<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.12.4 Is the procedure used to evaluate QC sample analyses in accordance with the most recent revision of Laboratory Modification LB-000029?		<input checked="" type="checkbox"/>	<input type="checkbox"/>	
7.13 Standard Operating Procedures (SOPs)				
7.13.1 Are the applicable laboratory SOPs available and followed by laboratory personnel (list)?		<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Document Title	Control No.	Description		
ASTM D5755-09	Rev. 3.1 (11/2/2012)			
TEM EDX	Rev. 1.2 (1/29/2014)			
EPA 100.2	Rev. 14.3(8/8/2013)			
ISO 10312	Rev. 9.1 (2/16/2012)			
ISO 13794	Rev. 1.2 (3/1/2013)			
SAED	Rev. 0 (12/7/2012)			
AHERA	Rev. 14.3 (4/5/2013)			
7.14 Document Control		Yes	No	Comments
7.14.1 Are all logbooks, notebooks, forms, or other laboratory documents legible, accurate, and complete (list)?		<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Document Title	Description/Comments			
Instrument Maintenance	TEM maintenance records			
Instrument Records	TEM calibration records			
Additional Comments:				

LIBBY-SPECIFIC ASBESTOS LABORATORY ON-SITE AUDIT CHECKLIST

Method: Not ApplicableDate(s) of On-site: March 25-26, 2014

8.0 POLARIZED LIGHT MICROSCOPY (PLM)		Yes	No	Comments
8.1	Are PLM areas adequate, clean, and orderly?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.2	Are steps taken to prevent the cross-contamination of equipment, supplies, and reagents?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Personnel Interviewed				
Name		Title		Experience
Garret Vliet		PLM Supervisor		4 Years
Melissa Klinedinst		PLM QC Group Leader		7 Years
8.3 Methods and Guidance Documents		Yes	No	Comments
8.3.1	Are the applicable guidance documents available for reference:			
8.3.1.1	EPA SOP SRC-Libby-01?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	All project specific documents and methods are available in the CDM Smith eRoom.
8.3.1.2	EPA SOP SRC-Libby-03?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.3.1.3	NIOSH 9002, Issue 2 - Asbestos (Bulk) by PLM?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.3.1.4	Others (list)?	<input type="checkbox"/>	<input checked="" type="checkbox"/>	
8.3.2	Are project-specific requirements communicated to laboratory personnel and available for reference:			
8.3.2.1	Laboratory Modifications:			
8.3.2.1.1	Most current revision of LB-000097?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.3.2.1.2	Most current revision of LB-000098?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.3.2.2	SAP Analytical Summaries?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.3.2.3	Project-specific Electronic Data Deliverables (EDDs)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.3.2.4	Others (list)?	<input type="checkbox"/>	<input checked="" type="checkbox"/>	
8.4 Equipment				
8.4.1	Ventilation Hoods:			
8.4.1.1	Checked routinely and recorded in a permanent logbook?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.4.2	Drying oven:			
8.4.2.1	Checked routinely and recorded in a permanent logbook?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	In bulk preparation area.
8.4.3	Muffle furnace:			
8.4.3.1	Checked routinely and recorded in a permanent logbook?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	In bulk preparation area.
Additional Comments:				

LIBBY-SPECIFIC ASBESTOS LABORATORY ON-SITE AUDIT CHECKLIST

Method: Not ApplicableDate(s) of On-site: March 25-26, 2014

8.0 POLARIZED LIGHT MICROSCOPY (PLM)	Yes	No	Comments
8.4.4 Analytical Balances:			
8.4.4.1 Do the balances meet the following criteria:			
8.4.4.1.1 Accurate to 1 mg (0.001 g)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.4.4.1.2 Upper range of at least 100 g?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.4.4.2 Checked routinely and recorded in a permanent logbook?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.4.4.2.1 Are balances calibrated using at least three weights?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.4.4.3 Has the balance been calibrated within the last 12 months by a third party vendor?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.4.4.4 Are the weights used traceable to national standards for weights and measures and certified by a third party with the last 5 years?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.4.4.5 Are both the balance and weights labeled with the following information:			Refer to Finding No. 4 of the On-site Audit Report.
8.4.4.5.1 Date of certification?	<input type="checkbox"/>	<input checked="" type="checkbox"/>	
8.4.4.5.2 Initials of individual performing the certification?	<input type="checkbox"/>	<input checked="" type="checkbox"/>	
8.4.4.5.3 Date next service is to be performed?	<input type="checkbox"/>	<input checked="" type="checkbox"/>	
8.5 Stereomicroscope			
8.5.1 Do stereomicroscopes meet the following requirements:			
8.5.1.1 Magnification range of 10X to 50X?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.5.1.2 Incandescent or fluorescent light source?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.6 Polarized Light Microscope			
8.6.1 Are PLMs equipped with the following:			
8.6.1.1 Light source and replacement bulbs?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.6.1.2 Binocular observation tube?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.6.1.3 Blue daylight filter?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.6.1.4 Oculars (10X)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.6.1.5 Objectives: 10X, 20X, and 40X (or similar)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.6.1.6 10X dispersion staining objective?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.6.1.7 A 360 degree graduated rotating stage?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.6.1.8 Polarizer and analyzer aligned at 90 degrees to one another?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.6.1.9 Bertrand lens?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.6.1.10 Substage condenser with iris diaphragm?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.6.1.11 Accessory slot for compensator plate?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.6.1.12 First order red (530-550 nanometer) compensator plate?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.6.1.13 Crosshair reticle?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.6.1.14 Adjustment tools?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Additional Comments:			

LIBBY-SPECIFIC ASBESTOS LABORATORY ON-SITE AUDIT CHECKLIST

Method: Not ApplicableDate(s) of On-site: March 25-26, 2014

8.0 POLARIZED LIGHT MICROSCOPY (PLM)			Yes	No	Comments
8.7 Are microscopes well-maintained, and are all routine and non-routine maintenance activities recorded in instrument-specific logbooks?			<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Instrument No.	Make	Model	Capabilities		
#2	Leica	DM750P	Standard		
#9	Leica	DMEP	Standard		
#10	Leica	DM750P	Standard		
8.8 Refractive Index Liquids			Yes	No	Comments
8.8.1 What refractive index liquids are available:					
8.8.1.1 High dispersion RI liquids from 1.620 to 1.640?			<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.8.1.2 1.550 high dispersion RI liquid?			<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.8.1.3 1.680 to 1.700 RI liquids?			<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.8.2 Are refractive index liquids checked daily for contamination?			<input checked="" type="checkbox"/>	<input type="checkbox"/>	Salt is used.
8.8.3 Are refractive index (RI) liquids calibrated monthly using a refractometer or other means (describe)?			<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.9 Reference Materials					
8.9.1 Does the laboratory maintain a library of Asbestos and non-Asbestos reference materials:					
8.9.1.1 NIST SRM 1866b (Ch, Am, and Cr)?			<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.9.1.2 NIST SRM 1867a (Tr, Ac, and An)?			<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.9.1.3 USGS LA PEs:			<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.9.1.3.1 LA 0.2% by mass?			<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.9.1.3.2 LA 1.0% by mass?			<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.9.1.3.3 Other (list)?			<input type="checkbox"/>	<input checked="" type="checkbox"/>	
8.9.1.4 Controlled LA Asbestos (USGS)?			<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.9.1.5 NIST testing round M12001 (Winchite/Richterite)?			<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.9.1.6 Non-Asbestos (i.e., Gypsum, Calcite, and Fiberglass)?			<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.10 PLM Calibration					
8.10.1 For PLM, is the following performed daily:					
8.10.1.1 Alignment?			<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.10.1.2 Stage and objectives centered?			<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.10.1.3 Optic axis centered?			<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.10.1.4 Alignment of the upper/lower polars?			<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.10.1.5 Centered through substage condenser and iris diaphragm?			<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.10.2 Microscope adjustments verified and recorded prior to sample analyses?			<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Additional Comments:					

LIBBY-SPECIFIC ASBESTOS LABORATORY ON-SITE AUDIT CHECKLIST

Method: Not ApplicableDate(s) of On-site: March 25-26, 2014

8.0 POLARIZED LIGHT MICROSCOPY (PLM)	Yes	No	Comments
8.11 PLM Analysis by NIOSH Method 9002:			
8.11.1 Does the laboratory perform PLM analyses on samples received from the Libby Superfund site? If "NO," proceed to Section 8.13 of the checklist.	<input type="checkbox"/>	<input checked="" type="checkbox"/>	
8.11.2 Are samples visually examined by stereomicroscope for the following: 8.11.2.1 Color? 8.11.2.2 Homogeneity? 8.11.2.3 Texture?	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	
8.11.3 Which of the following techniques are used to prepare samples for analysis: 8.11.3.1 Mortar & pestle? 8.11.3.2 Acid washing? 8.11.3.3 Ashing? 8.11.3.4 Solvents? 8.11.3.5 Other (list)?	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	
8.11.4 For non-friable, organically bound samples requiring ashing and/or acid reduction, are all necessary weights and tare weights measured and recorded?	<input type="checkbox"/>	<input type="checkbox"/>	
8.11.5 Are slides prepared using the appropriate refractive index liquid(s) and scanned for Asbestos fibers using the following optical properties: 8.11.5.1 Morphology? 8.11.5.2 Color? 8.11.5.3 Refractive indices? 8.11.5.4 Pleochroism? 8.11.5.5 Birefringence? 8.11.5.6 Extinction characteristics? 8.11.5.7 Sign of elongation? 8.11.5.8 Dispersion staining characteristics?	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	
8.11.6 Are the observed optical properties compared to Table 1 (Optical Properties of Asbestos Fibers) to determine the Asbestos mineral present?	<input type="checkbox"/>	<input type="checkbox"/>	
8.11.7 Is a qualitative assessment of Asbestos content made from both the gross and microscopic examinations?	<input type="checkbox"/>	<input type="checkbox"/>	
8.11.8 If no fibers are detected in a homogeneous samples are at least two additional slides prepared and analyzed prior to concluding no Asbestos is present?	<input type="checkbox"/>	<input type="checkbox"/>	
8.11.9 Is at least one optical property recorded for fibers determined to be non-Asbestos fibers?	<input type="checkbox"/>	<input type="checkbox"/>	
Additional Comments:			

LIBBY-SPECIFIC ASBESTOS LABORATORY ON-SITE AUDIT CHECKLIST

Method: Not ApplicableDate(s) of On-site: March 25-26, 2014

8.0 POLARIZED LIGHT MICROSCOPY (PLM)	Yes	No	Comments
8.13 PLM-VE (SOP SRC-Libby-03)			
8.13.1 Stereomicroscopic Examination:			
8.13.1.1 Are sample preparation activities performed within a HEPA-filtered hood?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.13.1.2 Is the entire sample transferred to an Asbestos-free container for examination?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.13.1.2.1 Is the container a minimum of 100 mm in diameter?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.13.1.2.2 If non-disposable, is a daily contamination check performed?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.13.1.3 Is the entire sample examined for homogeneity and the presence of suspect fibers?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.13.1.4 Are suspect fibers removed with fine forceps and mounted in the appropriate RI liquid for PLM analysis?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.13.1.5 Are the following stereomicroscopic findings recorded:			
8.13.1.5.1 Sample appearance?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.13.1.5.2 Estimated percentage of LA?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.13.1.5.3 Estimated percentage of other Asbestos types?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.13.2 Determination for Ashing the Sample:			
8.13.2.1 Are soil samples containing a significant amount of artifacts ashed prior to being prepared for random PLM mounts?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.13.2.1.1 Are samples ashed in a muffle furnace at approximately 480 °C?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.13.2.1.2 Are the necessary gravimetric measurements recorded for the determination of "Pre-ash percent Asbestos"?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.13.3 Determination for Additional Grinding:			
8.13.3.1 As necessary, are samples ground by a mortar and pestle?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.13.3.2 Is the mortar and pestle cleaned between samples and a daily contamination check performed?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.13.4 Slide Preparation for PLM-VE:			
8.13.4.1 Are a minimum of five random sub-samples mounted in the appropriate RI liquid (1.620-1.640) for measurement of LA optical properties?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.13.5 Supplemental Stereomicroscopic Evaluation:			
8.13.5.1 Following the random slide mount preparation, is the container agitated to cause the particulate to settle and Asbestos fibers sort to the surface?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.13.5.2 Is the sample re-examined, and the fiber pick procedure repeated?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Additional Comments: Analyst demonstrated modified (LB-000096) PLM-VE technique on previously analyzed Inter-lab samples that were weakly discordant from original results. For each of the samples reanalyzed, the result reported for the inter-laboratory study was confirmed.			

LIBBY-SPECIFIC ASBESTOS LABORATORY ON-SITE AUDIT CHECKLIST

Method: Not ApplicableDate(s) of On-site: March 25-26, 2014

8.0 POLARIZED LIGHT MICROSCOPY (PLM)	Yes	No	Comments
8.14 PLM-GRAV (SOP SRC-Libby-01)			
8.14.1 Stereomicroscopic Examination:			
8.14.1.1 Is the entire sample weighed and placed in an appropriate container?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.1.2 Does the stereomicroscopic examination include:			
8.14.1.2.1 Examination of multiple fields of view over the entire sample?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.1.2.2 Probing of the sample and breaking clumps where possible?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.1.2.3 Manipulation of the sample with the appropriate tools?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.1.2.4 Observations for homogeneity, texture, friability, color, and extent of any Asbestos content?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.1.3 Does the analyst refrain from segregating and weighing particles smaller than 2-3 mm (1/10 inch)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.1.4 If no particles larger than 2-3 mm or larger are present, are one of the following recorded:			
8.14.1.4.1 No Asbestos detected (ND)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.1.4.2 Trace levels of Asbestos observed, but not quantified (Tr)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.2 Examination by PLM:			
8.14.2.1 Are tentatively identified Asbestos particles examined by PLM as described in SOP SRC-Libby-3 (Section 8.12 of this checklist)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.2.2 If Asbestos particles are determined to be OA, are they further characterized:			
8.14.2.2.1 Amosite (AMOS)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.2.2.2 Anthophyllite (ANTH)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.2.2.3 Crocidolite (CROC)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.2.2.4 Unknown (UNK)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.2.3 Is the total weight of each type of positively identified Asbestos measured and recorded?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.3 Record Keeping:			
8.14.3.1 Is the data log sheet provided in Attachment 1 of the SOP used to record weights the initial (coarse fraction) and segregated Asbestos?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Additional Comments:			

LIBBY-SPECIFIC ASBESTOS LABORATORY ON-SITE AUDIT CHECKLIST

Method: Not ApplicableDate(s) of On-site: March 25-26, 2014

8.0 POLARIZED LIGHT MICROSCOPY (PLM)	Yes	No	Comments
8.14.4 Classification of Asbestos Mineral Type:			
8.14.4.1 Using PLM, is entire area of each prepared slide examined for Asbestos, non-Asbestos and matrix material?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.4.2 Is positive identification determined from the following six optical properties:			
8.14.4.2.1 Habit (Asbestos or non-Asbestos)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.4.2.2 Color & Pleochroism (if present)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.4.2.3 Both Alpha and Gamma Refractive indices?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.4.2.4 Birefringence?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.4.2.5 Extinction angle?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.4.2.6 Sign of elongation (positive-slow or negative fast)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.4.3 Based on the optical properties, is Asbestos classified into one of the following three categories:			
8.14.4.3.1 Libby Amphibole (LA)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.4.3.2 Other Amphibole (OA)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.4.3.3 Chrysotile (CH)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.4.4 Is at least one optical property recorded for observed non-Asbestos fibers?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.5 Quantification of Asbestos Content:			
8.14.5.1 Is Asbestos reported as either mass or area percent for LA?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.5.2 Are other, non-LA, Asbestos types reported in area percent?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.5.3 Are reference materials used to aid in visual estimation:			
8.14.5.3.1 LA PE reference materials (0.2% or 1.0%)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.5.3.2 Are visual estimates of greater than 1% LA performed using calibration standards made in-house from NIST SRMs and NIST PEs?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.6 Are calibrated visual estimates determined from both the detailed stereomicroscopic observations and examination of the total area for all five random slide mounts?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.7 Are LA results reported in the appropriate bin categories (PLM-VE only):			
8.14.7.1 Non-detects recorded as Bin A?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.7.2 Less than 0.2% LA recorded as Bin B1?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.7.3 Greater than 0.2%, but less than 1% LA recorded as Bin B2?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.14.7.4 Equal to or greater than 1% LA recorded as Bin C, with the percentage recorded as a whole number?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Additional Comments:			

LIBBY-SPECIFIC ASBESTOS LABORATORY ON-SITE AUDIT CHECKLIST

Method: Not ApplicableDate(s) of On-site: March 25-26, 2014

8.0 POLARIZED LIGHT MICROSCOPY (PLM)		Yes	No	Comments
8.15 Quality Control Analyses				
8.15.1 Are the following types of QC analyses performed at the required frequencies:				
8.15.1.1 Laboratory Duplicate Self-check (LDS) at a frequency of 2%?	<input checked="" type="checkbox"/>	<input type="checkbox"/>		
8.15.1.2 Laboratory Duplicate Cross-check (LDC) at the correct frequency:				
8.15.1.2.1 At a frequency of 4% for PLM-VE?	<input checked="" type="checkbox"/>	<input type="checkbox"/>		
8.15.1.2.2 At a frequency of 8% for PLM-Grav?	<input checked="" type="checkbox"/>	<input type="checkbox"/>		
8.15.1.3 Laboratory Duplicate Cross-check re-preparation (LDCR) at a frequency of 4% (PLM-VE only)?	NA	NA		This requirement has not been made final yet.
8.15.1.4 Is a second analyst available for LDC and LDCR analyses?	<input checked="" type="checkbox"/>	<input type="checkbox"/>		
8.15.2 For sample containing LA, are LDS, LDC and LDCR analyses considered acceptable if:				
8.15.2.1 LA results are within 1 Bin category?	<input checked="" type="checkbox"/>	<input type="checkbox"/>		
8.15.2.2 LA results are ≤1% LA?	<input checked="" type="checkbox"/>	<input type="checkbox"/>		
Note: For LA results greater than 1%, the laboratory should refer to their internal QA/QC system.				
8.15.3 Is the appropriate corrective action taken when LDS, LDC, or LDCR analyses do not meet acceptance criteria (describe)?		<input checked="" type="checkbox"/>	<input type="checkbox"/>	
8.16 Standard Operating Procedures (SOPs)				
8.16.1 Are the applicable laboratory SOPs available and followed by laboratory personnel (list)?		<input type="checkbox"/>	<input checked="" type="checkbox"/>	Refer to Finding No. 4 of the On-site Audit Report.
Document Title	Control No.	Description		
PLM SOP	Rev. 12.1 (10/15/2012)			
RI Liquid Calibration	Rev. 1.4 (9/2/2010)			
8.17 Document Control		Yes	No	Comments
8.17.1 Are all logbooks, notebooks, forms, or other laboratory documents legible, accurate, and complete (list)?		<input checked="" type="checkbox"/>	<input type="checkbox"/>	
Document Title	Description/Comments			
Balance Calibration Logbook				
Scope Calibration/Contamination Logbook	For Scope #10			
Additional Comments:				

LIBBY-SPECIFIC ASBESTOS LABORATORY ON-SITE AUDIT CHECKLIST

Method: Not ApplicableDate(s) of On-site: March 25-26, 2014

9.0 DATA MANAGEMENT	PCM	TEM	PLM	Comments
9.1 Data Package Review and Assembly				
9.1.1 Are deliverables reviewed to ensure project-specific requirements in the following are met:				
9.1.1.1 Request for Modifications to Laboratory Activities?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	All project specific documents and EDD templates are available in the CDM Smith eRoom.
9.1.1.2 Project-specific SOPs?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
9.1.1.3 SAP Analytical Summaries?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
9.1.1.4 Project-specific Electronic Data Deliverables (EDDs)?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
9.1.1.5 Other (list)? _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
9.1.2 Do scanned data packages meet the following requirements:				
9.1.2.1 A Case Narrative briefly describing the analytical method(s), any deviations from the method, and any other discrepancies?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
9.1.2.2 A copy of the signed COC and other documentation included with the COC?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
9.1.2.3 Sample information and final analytical results for all sample and QC analyses?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
9.1.2.4 Copies of handwritten bench sheets and sample preparation logs for indirect preparations?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
9.1.2.5 TEM spectra?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
9.1.2.6 TEM structure morphology photos (if required)?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
9.1.2.7 Additional documentation relevant to the test report (e.g., email)?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
9.1.2.8 Reports are paginated?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
9.1.3 Are all deliverables reviewed for completeness and accuracy prior to being submitted:				
9.1.3.1 Hardcopy deliverables?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
9.1.3.2 Electronic deliverables?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
9.1.4 Are all reviews documented?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
9.1.5 Are supplements and corrections to minor errors (that do not affect the results or EDD) submitted in the following manner:				
9.1.5.1 A description of the corrections made (can be provided on the cover page or in a case narrative)?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
9.1.5.2 The name, title and signature of an approved signatory?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
9.1.5.3 Any additional paperwork showing corrections?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
9.1.5.4 Verification initials/date and validation initials/date?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
9.1.5.5 Update page numbers to include the added supplement?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
9.1.6 A new, complete test report is generated anytime data are affected (resulting in a change to the EDD) and new data loaded into Scribe?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
Additional Comments				

LIBBY-SPECIFIC ASBESTOS LABORATORY ON-SITE AUDIT CHECKLIST

Method: Not ApplicableDate(s) of On-site: March 25-26, 2014

9.0 DATA MANAGEMENT	PCM	TEM	PLM	Comments
9.2 Data Submission				
9.2.1 Is the date that the electronic data deliverables (EDDs) are uploaded to the FTP site tracked and recorded?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
9.2.2 Are EDDs also e-mailed to the required recipients?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
9.2.3 Is the date that scanned data deliverables are uploaded to the FTP site tracked and recorded?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
9.2.4 Is a system in place to ensure deliverables are submitted within the required TATS:				A spreadsheet is utilized.
9.2.4.1 Electronic Deliverables?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	Refer to Finding No. 5 of the On-site Audit Report.
9.2.4.2 Hardcopy (scanned)?	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
9.2.5 Does the laboratory post the required PLM, TEM, and PCM calibration data to the CDM eRoom as specified in the most recent revision of Laboratory Modification LB-000085?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
9.3 Data Storage and Archiving				
9.3.1 Are electronic files archived onto suitable media on a frequent basis? How often? <u>Daily</u>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	Refer to Section 5.10.2 of the Control of Records SOP (Rev. 3.1, 2/04/2014).
9.3.2 Are all hardcopy data stored in a secured location with limited access (e.g., locking file cabinet, etc.)?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
Additional Comments: Note that, for the Libby Superfund project, all data entry, data review, and data package assembly procedures are performed for all of the EMSL Analytical laboratories in Cinnaminson, NJ by the Special Projects Data Coordinator.				

LIBBY-SPECIFIC ASBESTOS LABORATORY ON-SITE AUDIT CHECKLIST

Method: Not ApplicableDate(s) of On-site: March 25-26, 2014

10.0 QUALITY ASSURANCE/QUALITY CONTROL	PCM	TEM	PLM	Comments
10.1 Laboratory Certifications				
10.1.1 Is the laboratory accredited for Asbestos analysis under the National Voluntary Laboratory Accreditation Program (NVLAP):				Laboratory ID: 101048-0 Issued: 07/01/2013 Expire: 06/30/2014
10.1.1.1 Asbestos Fiber Analysis (TEM Method)?	NA	<input checked="" type="checkbox"/>	NA	
10.1.1.2 Asbestos Fiber Analysis (PLM Method)?	NA	NA	<input checked="" type="checkbox"/>	
10.1.2 Is the laboratory accredited for Asbestos analysis under the American Industrial Hygiene Association (AIHA), and does it participate in the National Institute for Occupational Safety and Health (NIOSH) Proficiency Analytical Testing (PAT) Program?	<input checked="" type="checkbox"/>	NA	NA	Laboratory ID: 100194 Issued: 07/31/2012 Expire: 07/01/2013
10.2 Training				
10.2.1 Have all analysts undergone training on the proper usage of the equipment and instrumentation used in the respective areas?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	Refer to Finding No. 3 of the On-site Audit Report.
10.2.2 Have all analysts demonstrated proficiency through the preparation and/or analysis of standards or samples of known values?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
10.2.3 Are training records maintained in analyst-specific files?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
10.3 Internal Audits				
10.3.1 Are internal audits conducted on an annual basis using an appropriate checklist?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	Conducted 05/08-10/2013.
10.3.1.1 Are internal audit reports available for review?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	Refer to Finding No. 6 of the On-site Audit Report.
10.4 Corrective/Preventive Action:				
10.4.1 Can the laboratory demonstrate the sequence of problem identification, corrective action, and resumption of duties?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
10.5 Quality Records				
10.5.1 Are SOPs available in the applicable areas for all laboratory-specific procedures?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
10.5.2 Does the laboratory have a Quality Assurance Manual?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
10.5.3 Does the laboratory compile monthly quality assurance/quality control reports?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
10.6 Environmental Controls/Laboratory Monitoring				
10.6.1 Does the laboratory conduct an environmental monitoring program?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
10.6.2 Is quarterly air monitoring performed in all laboratory areas?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
10.6.2.1 Are the collected samples analyzed by TEM with a target analytical sensitivity of 0.005 structures/cc?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	Air monitoring results for 2013/2014 reviewed prior to audit, with no structures detected.
10.6.2.2 If LA is detected, are the affected areas thoroughly cleaned and a new set of samples collected and analyzed?	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
Additional Comments:				