



UNITED STATES ENVIRONMENTAL PROTECTION AGENCY  
REGION 5 CHICAGO REGIONAL LABORATORY  
536 SOUTH CLARK STREET (ML-10C)  
CHICAGO, ILLINOIS 60605

**ELECTRONIC LABORATORY DATA PACKAGE  
FINAL LEVEL IV  
COVER PAGE**

LIMS Work order(s):	
Analysis:	
Primary Analyst:	
Date:	
Data Reporting Qualtrax Workflow ID:	

Digital Signature of Primary Analyst:

**Digital Signature Agreement:** By signing above the primary analyst understands and agrees that they will be held legally bound, obligated, and responsible for the use of their digital signature as they would be by using their hand-written signature.



UNITED STATES ENVIRONMENTAL PROTECTION AGENCY  
REGION 5 CHICAGO REGIONAL LABORATORY  
536 SOUTH CLARK STREET  
CHICAGO, ILLINOIS 60605

**Date:** 7/19/2018  
**Subject:** Review of Region 5 Data for Milliken Chemical Dewey Plant  
**To:** Region 4  
980 College Station Road  
Athens, GA 30605  
**From:** Colin Breslin, Chemist  
US EPA Region 5 Chicago Regional Laboratory

The data transmitted under this cover memo successfully passed CRL's data review procedures as documented in the current Quality Management Plan and applicable Standard Operating Procedures. In accordance with the EPA QA/G-8 *Guidance on Environmental Data Verification and Data Validation* and the U.S. EPA Region 5 RMD QMP, CRL performs data verification on all the data generated internally. CRL does not perform data validation or quality assessment procedures.

This report was reviewed and the information provided herein accurately represents the analysis performed.

A handwritten signature in black ink that reads "Colin Breslin". The signature is written in a cursive style and is positioned above a horizontal line.

*Please contact the analyst with any technical report issues, Robert Thompson at (312)-353-9078 for sample project concerns, and Sylvia Griffin at (312)-353-9073 with data transmittal questions. Thank you.*

**Attached are Results for: Milliken Chemical Dewey Plant**

**Analyses included in this report:**

Ignitability by Setaflash

Water content, Karl Fisher Titration

Report Name: 1806004 Ignitability by Setaflash Water content, Karl Fisher Titration FINAL Jul 19 18 1448



Environmental Protection Agency Region 5  
**Chicago Regional Laboratory**

536 South Clark Street, Chicago, IL 60605  
Phone:(312)353-8370 Fax:(312)886-2591

Region 4  
980 College Station Road  
Athens GA, 30605

Project: Milliken Chemical Dewey Plant  
Project Number: [none]  
Project Manager: Jeffrey Hendel

**Reported:**  
Jul-19-18 14:48

**Accredited Analyses included in this Report**



**Method:** *EPA 1020B in Water*  
**Analysis:** *Ignitability by Setaflash*  
**Analyte**

**Certifications**

Ignitability by Flashpoint ISO/IEC 17025:2005

**Method:** *EPA 9000 in Water*  
**Analysis:** *Water content,Karl Fisher Titration*  
**Analyte**

**Certifications**

Water content ISO/IEC 17025:2005

Analytes not listed above are not accredited by ANAB.



Environmental Protection Agency Region 5  
**Chicago Regional Laboratory**

536 South Clark Street, Chicago, IL 60605  
Phone:(312)353-8370 Fax:(312)886-2591

Region 4  
980 College Station Road  
Athens GA, 30605

Project: Milliken Chemical Dewey Plant  
Project Number: [none]  
Project Manager: Jeffrey Hendel

**Reported:**  
Jul-19-18 14:48

## Analysis Case Narrative

### General Information

Two samples for the analysis of Ignitability by Setaflash (flashpoint) and Water Content, Karl Fisher Titration were received at the Chicago Regional Laboratory (CRL) on June 06, 2018. Holding times do not apply to these samples. The designated analyst, Colin Breslin, can be reached at 312-886-2912.

The flashpoint samples were prepared and analyzed according to CRL SOP AIG048A Version 5 (based on SW-846: 1020B).

The water content samples were prepared and analyzed according to CRL SOP AIG015A Version 4 (based on SW-846: 9000).

### Sample Analysis and Results

The data reported herein meets the requirements of the CRL analytical SOP used for analysis. There were no laboratory specifications stated in the Milliken Chemical Dewey Plant Quality Assurance Project Plan (QAPP), dated September 19, 2017. The QAPP is saved on the CRL share drive with file name: Miliken\_QAPP\_17-0496.pdf.

Water content results were reported as not detected at the CRL analytical SOP reporting limit of 0.5% water content, which showed the samples were essentially non-aqueous.

### Quality Control

All Quality Control (QC) audits were within CRL limits for the requested analytes or did not result in qualification of the data.



Environmental Protection Agency Region 5  
**Chicago Regional Laboratory**

536 South Clark Street, Chicago, IL 60605  
Phone:(312)353-8370 Fax:(312)886-2591

**WORK ORDER**

Printed: 6/7/2018 11:07:31AM

**1806004**

**US EPA Region 5 Chicago Regional Laboratory**

**Client:** Region 4  
**Project:** Milliken Chemical Dewey Plant

**Project Manager:** Angela Ockrassa Davis  
**Project Number:** [none]

**Report To:**

Jeffrey Hendel  
Region 4

980 College Station Road  
Athens, GA 30605

Phone: (706) 355-8839  
Fax:

Date Due: Jul-23-18 15:00 (46 day TAT)

Received By: Robert Snyder

Date Received: Jun-06-18 09:15

Logged In By: Robert Snyder

Date Logged In: Jun-06-18 13:57

Samples Received at: 22.8 °C  
Sample tags/labels Yes  
Seals Intact Yes  
Received on ice No  
Paperwork Included Yes

**Work Order Comments:**  
Copy/Relog from 1805008.

**Sample ID:** [1806004-01](#) **Sampled:** [Sep-26-17 12:25](#) **Matrix:** [Water](#)

**Sample Name:** [01 Tetramer](#) **Sample Location/Comments:** [E173905-01](#)

**Sample Comments:**

Analysis	Hold time (days)	Expires	Comments
Ignitability by Setaflash	365	Sep-26-18 12:25	pH = 4
Water content,Karl Fisher Titration	365	Sep-26-18 12:25	pH = 4

**Sample ID:** [1806004-02](#) **Sampled:** [Sep-26-17 12:35](#) **Matrix:** [Water](#)

**Sample Name:** [02 Tetramer](#) **Sample Location/Comments:** [E173905-02](#)

**Sample Comments:**

Analysis	Hold time (days)	Expires	Comments
Ignitability by Setaflash	365	Sep-26-18 12:35	pH = 4
Water content,Karl Fisher Titration	365	Sep-26-18 12:35	pH = 4

WORK ORDER MEMO: Samples rec'd at ambient temperature. Last sample set for this project was also, and NCR 14124 was opened. NCR was resolved as not being an issue, as flash point analysis does not have a requirement to be cooled but the sample requirements table just hadn't been updated yet. Update currently pending, so no NCR will be opened for this shipment.

**REVIEWED**

By Robert Thompson at 11:08 am, Jun 07, 2018

Reviewed By

Date



Environmental Protection Agency Region 5  
**Chicago Regional Laboratory**

536 South Clark Street, Chicago, IL 60605  
Phone:(312)353-8370 Fax:(312)886-2591

Region 4  
980 College Station Road  
Athens GA, 30605

Project: Milliken Chemical Dewey Plant  
Project Number: [none]  
Project Manager: Jeffrey Hendel

**Reported:**  
Jul-19-18 14:48

**Water Content by SW846 Method 9000**  
**US EPA Region 5 Chicago Regional Laboratory**

**01 Tetramer (1806004-01) Matrix: Water Sampled: Sep-26-17 12:25 Received: Jun-06-18 09:15**

Analyte	Result	Flags / Qualifiers	MDL	Reporting Limit	Units	Dilution	Batch	Prepared	Analyzed
Water content	U			0.5	%	1	B18F011	Jun-12-18	Jun-12-18

**02 Tetramer (1806004-02) Matrix: Water Sampled: Sep-26-17 12:35 Received: Jun-06-18 09:15**

Analyte	Result	Flags / Qualifiers	MDL	Reporting Limit	Units	Dilution	Batch	Prepared	Analyzed
Water content	U			0.5	%	1	B18F011	Jun-12-18	Jun-12-18



Environmental Protection Agency Region 5  
**Chicago Regional Laboratory**

536 South Clark Street, Chicago, IL 60605  
Phone:(312)353-8370 Fax:(312)886-2591

Region 4  
980 College Station Road  
Athens GA, 30605

Project: Milliken Chemical Dewey Plant  
Project Number: [none]  
Project Manager: Jeffrey Hendel

**Reported:**  
Jul-19-18 14:48

**Ignitability, Flash Point, EPA 1020B (modified)**  
**US EPA Region 5 Chicago Regional Laboratory**

**01 Tetramer (1806004-01)**

**Matrix: Water**

**Sampled: Sep-26-17 12:25**

**Received: Jun-06-18 09:15**

Analyte	Result	Flags / Qualifiers	MDL	Reporting Limit	Units	Dilution	Batch	Prepared	Analyzed
<b>Ignitability by Flashpoint</b>	<b>136</b>				Degrees F	1	B18F004	Jun-08-18	Jun-08-18

**02 Tetramer (1806004-02)**

**Matrix: Water**

**Sampled: Sep-26-17 12:35**

**Received: Jun-06-18 09:15**

Analyte	Result	Flags / Qualifiers	MDL	Reporting Limit	Units	Dilution	Batch	Prepared	Analyzed
<b>Ignitability by Flashpoint</b>	<b>136</b>				Degrees F	1	B18F004	Jun-08-18	Jun-08-18



Environmental Protection Agency Region 5  
**Chicago Regional Laboratory**

536 South Clark Street, Chicago, IL 60605  
Phone:(312)353-8370 Fax:(312)886-2591

Region 4  
980 College Station Road  
Athens GA, 30605

Project: Milliken Chemical Dewey Plant  
Project Number: [none]  
Project Manager: Jeffrey Hendel

**Reported:**  
Jul-19-18 14:48

**Notes and Definitions**

U Not Detected  
NR Not Reported  
Q QC limit Exceeded





Environmental Protection Agency Region 5  
**Chicago Regional Laboratory**

536 South Clark Street, Chicago, IL 60605  
Phone:(312)353-8370 Fax:(312)886-2591

## Quality Control Summary

Date: 7/19/2018

Analyst: Colin Breslin

Project: Milliken Chemical Dewey Plant

Analyte: Ignitability by Flashpoint

Analysis: Ignitability by Setaflash

QC Sample Name	Result	Units	Flags / Qualifiers	MDL	Reporting Limit	Spike Level	Source Result	%REC	%REC Limits	% RPD	% RPD Limit
B18F004-DUP2 Source: 1806004-01	136	Degrees F					136			0.0147	2.23
B18F004-SRM1	80.0	Degrees F				81.00		98.8%	98.1-101.9		
B18F004-SRM2	81.0	Degrees F				81.00		100%	98.1-101.9		
B18F004-SRM3	80.6	Degrees F	CustomFlag			81.00		99.4%	98.1-101.9		
B18F004-SRM4	81.0	Degrees F				81.00		100%	98.1-101.9		
B18F004-SRM5	81.0	Degrees F				81.00		100%	98.1-101.9		
B18F004-SRM6	81.0	Degrees F	CustomFlag			81.00		100%	98.1-101.9		
B18F004-SRM7	80.8	Degrees F				81.00		99.7%	98.1-101.9		
B18F004-SRM8	79.7	Degrees F				81.00		98.4%	98.1-101.9		
B18F004-SRM9	80.2	Degrees F	CustomFlag			81.00		99.1%	98.1-101.9		

Analyte: Water content

Analysis: Water content,Karl Fisher Titration

QC Sample Name	Result	Units	Flags / Qualifiers	MDL	Reporting Limit	Spike Level	Source Result	%REC	%REC Limits	% RPD	% RPD Limit
B18F011-DUP7 Source: 1806004-01	0.01	%			0.5		0.01			0.00	10
B18F011-DUP8 Source: 1806004-02	0.01	%			0.5		0.01			0.00	10
B18F011-SRM1	1.04	%				1.000		104%	97-107		
B18F011-SRM2	1.04	%				1.000		104%	97-107		
B18F011-SRM3	1.03	%				1.000		103%	97-107		
B18F011-SRM4	1.04	%				1.000		104%	97-107		



Environmental Protection Agency Region 5  
**Chicago Regional Laboratory**

536 South Clark Street, Chicago, IL 60605  
Phone:(312)353-8370 Fax:(312)886-2591

---

## **Quality Control Summary**

---

**Date:** 7/19/2018

**Analyst:** Colin Breslin

**Project:** Milliken Chemical Dewey Plant

---

### **Notes and Definitions**

CustomFlag B18F004-SRM9: average of B18F004-SRM7 & B18F004-SRM8

CustomFlag B18F004-SRM6: average of B18F004-SRM4 & B18F004-SRM5

CustomFlag B18F004-SRM3: average of B18F004-SRM1 & B18F004-SRM2

U Not Detected

NR Not Reported

Q QC limit Exceeded

***REVIEWED***

***By Colin Breslin at 2:44 pm, Jul 19, 2018***

---

Colin Breslin, Chemist



Environmental Protection Agency Region 5  
**Chicago Regional Laboratory**

536 South Clark Street, Chicago, IL 60605  
Phone:(312)353-8370 Fax:(312)886-2591

**WORK ORDER**

Printed: 6/7/2018 11:07:31AM

**1806004**

**US EPA Region 5 Chicago Regional Laboratory**

**Client:** Region 4 **Project Manager:** Angela Ockrassa Davis  
**Project:** Milliken Chemical Dewey Plant **Project Number:** [none]

**Report To:**

Jeffrey Hendel  
Region 4

980 College Station Road  
Athens, GA 30605

Phone: (706) 355-8839  
Fax:

Date Due: Jul-23-18 15:00 (46 day TAT)

Received By: Robert Snyder

Date Received: Jun-06-18 09:15

Logged In By: Robert Snyder

Date Logged In: Jun-06-18 13:57

Samples Received at: 22.8 °C  
Sample tags/labels Yes  
Seals Intact Yes  
Received on ice No  
Paperwork Included Yes

**Work Order Comments:**  
Copy/Relog from 1805008.

**Sample ID:** [1806004-01](#) **Sampled:** [Sep-26-17 12:25](#) **Matrix:** [Water](#)

**Sample Name:** [01 Tetramer](#) **Sample Location/Comments:** [E173905-01](#)

**Sample Comments:**

Analysis	Hold time (days)	Expires	Comments
Ignitability by Setaflash	365	Sep-26-18 12:25	pH = 4
Water content,Karl Fisher Titration	365	Sep-26-18 12:25	pH = 4

**Sample ID:** [1806004-02](#) **Sampled:** [Sep-26-17 12:35](#) **Matrix:** [Water](#)

**Sample Name:** [02 Tetramer](#) **Sample Location/Comments:** [E173905-02](#)

**Sample Comments:**

Analysis	Hold time (days)	Expires	Comments
Ignitability by Setaflash	365	Sep-26-18 12:35	pH = 4
Water content,Karl Fisher Titration	365	Sep-26-18 12:35	pH = 4

WORK ORDER MEMO: Samples rec'd at ambient temperature. Last sample set for this project was also, and NCR 14124 was opened. NCR was resolved as not being an issue, as flash point analysis does not have a requirement to be cooled but the sample requirements table just hadn't been updated yet. Update currently pending, so no NCR will be opened for this shipment.

**REVIEWED**

By Robert Thompson at 11:08 am, Jun 07, 2018

Reviewed By

Date

## Breslin, Colin

---

**From:** Qualtrax <noreply@qualtraxcloud.com>  
**Sent:** Tuesday, July 17, 2018 8:19 AM  
**To:** Breslin, Colin  
**Subject:** Notification of Client Contact "WO 1806004 Flashpoint" (ID 15093)

The following is a summary of Client Contact Workflow ID 15093. Adjust your project plan accordingly.

Contact Date	Client Name and Affiliation	Project and/or Work Order	Contact Description	Action Request
5/30/2018 12:00:00 AM	US EPA R4 / Jeffrey Hendel	Milliken Chemical Dewey Plant / 1806004	request to send flashpoint samples.	agreed to accept samples.

## Breslin, Colin

---

**From:** Awanya, Francis  
**Sent:** Wednesday, May 30, 2018 2:05 PM  
**To:** Thompson, Robert A.; Breslin, Colin; Persoon, Carolyn; Schupp, George  
**Subject:** RE: request form R4 flashpoint

**Follow Up Flag:** Follow up  
**Flag Status:** Flagged

Okay Rob.

Francis A. Awanya  
Chemist  
Chicago Regional Laboratory  
USEPA Region 5  
536 South Clark Street (ML-10C)  
Chicago, IL 60605  
Telephone: (312)886-3682  
Fax: (312)692-2404

---

**From:** Thompson, Robert A.  
**Sent:** Wednesday, May 30, 2018 1:55 PM  
**To:** Breslin, Colin <breslin.colin@epa.gov>; Awanya, Francis <awanya.francis@epa.gov>; Persoon, Carolyn <persoon.carolyn@epa.gov>; Schupp, George <schupp.george@epa.gov>  
**Subject:** FW: request form R4 flashpoint

2 flashpoint samples from R4 requested. They wanted 30-days. Okay to accept? 250mL poly?



Robert Thompson, MPH  
Chemist | Sample Coordinator  
U.S.EPA R5 Chicago Regional Laboratory  
536 S. Clark St. ML-10C | Chicago, IL 60605  
(312) 353-9078 (Direct) | (312) 353-9096 (Fax)



Think Green - Not every email needs to be printed

---

**From:** Hendel, Jeffrey  
**Sent:** Wednesday, May 30, 2018 1:49 PM  
**To:** Thompson, Robert A. <[Thompson.Robert@epa.gov](mailto:Thompson.Robert@epa.gov)>  
**Subject:** RE: request form

Here is the completed form. Please let me know how much sample that the lab requires for performing flashpoint.

Thanks for your help,  
Jeff

Jeffrey R. Hendel

Chief, Inorganic Chemistry Section  
U.S. Environmental Protection Agency  
980 College Station Road  
Athens, GA 30605  
Direct: 706-355-8839  
Cell: 706-340-2145  
[hendel.jeffrey@epa.gov](mailto:hendel.jeffrey@epa.gov)

---

**From:** Thompson, Robert A.  
**Sent:** Wednesday, May 30, 2018 2:08 PM  
**To:** Hendel, Jeffrey <[Hendel.Jeffrey@epa.gov](mailto:Hendel.Jeffrey@epa.gov)>  
**Subject:** request form



Robert Thompson, MPH  
Chemist | Sample Coordinator  
U.S.EPA R5 Chicago Regional Laboratory  
536 S. Clark St. ML-10C | Chicago, IL 60605  
(312) 353-9078 (Direct) | (312) 353-9096 (Fax)




























Think Green - Not every email needs to be printed

**Work Order: 1806004****Analysis: Water Content, Karl Fisher Titration**

I:\Projects (New File Structure)\Milliken Chemical Dewey Plant\1.

A&amp;I\CBreslin\1806004\KFT\Water Content

Milliken Chemical Dewey Plant > 1. A&I > CBreslin > 1806004 > KFT > Water Content					Search Water Content	
Name		Date modified	Type	Size		
	1805008-01_D0505ACRLKTF_20180612-122656.mdet	6/12/2018 12:26 PM	MDET File	77 KB		
	1805008-01_D0505ACRLKTF_20180612-122657.mdet	6/12/2018 12:26 PM	MDET File	77 KB		
	1805008-02_D0505ACRLKTF_20180612-122658.mdet	6/12/2018 12:27 PM	MDET File	78 KB		
	1805008-03_D0505ACRLKTF_20180612-122656.mdet	6/12/2018 12:26 PM	MDET File	78 KB		
	1805008-03_D0505ACRLKTF_20180612-122656_1.mdet	6/12/2018 12:26 PM	MDET File	78 KB		
	1805008-04_D0505ACRLKTF_20180612-122655.mdet	6/12/2018 12:26 PM	MDET File	78 KB		
	1805008-05_D0505ACRLKTF_20180612-122655.mdet	6/12/2018 12:26 PM	MDET File	78 KB		
	1805008-05_D0505ACRLKTF_20180612-122655_1.mdet	6/12/2018 12:26 PM	MDET File	78 KB		
	1805008-06_D0505ACRLKTF_20180612-122654.mdet	6/12/2018 12:26 PM	MDET File	78 KB		
	1806004-01_D0505ACRLKTF_20180612-122658.mdet	6/12/2018 12:27 PM	MDET File	77 KB		
	1806004-01_D0505ACRLKTF_20180612-122658_1.mdet	6/12/2018 12:27 PM	MDET File	77 KB		
	1806004-02_D0505ACRLKTF_20180612-122657.mdet	6/12/2018 12:26 PM	MDET File	77 KB		
	B18F011-DUP1_D0505ACRLKTF_20180612-122656.mdet	6/12/2018 12:26 PM	MDET File	77 KB		
	B18F011-DUP2_D0505ACRLKTF_20180612-122658.mdet	6/12/2018 12:27 PM	MDET File	77 KB		
	B18F011-DUP3_D0505ACRLKTF_20180612-122655.mdet	6/12/2018 12:26 PM	MDET File	77 KB		
	B18F011-DUP4_D0505ACRLKTF_20180612-122655.mdet	6/12/2018 12:26 PM	MDET File	77 KB		
	B18F011-DUP5_D0505ACRLKTF_20180612-122654.mdet	6/12/2018 12:26 PM	MDET File	77 KB		
	B18F011-DUP6_D0505ACRLKTF_20180612-122654.mdet	6/12/2018 12:26 PM	MDET File	77 KB		
	B18F011-DUP7_D0505ACRLKTF_20180612-122657.mdet	6/12/2018 12:27 PM	MDET File	77 KB		
	B18F011-DUP8_D0505ACRLKTF_20180612-122657.mdet	6/12/2018 12:26 PM	MDET File	77 KB		
	Hydranal 10.0 (B18F011-SRM1)_D0505ACRLKTF_20180612-122659.mdet	6/12/2018 12:27 PM	MDET File	78 KB		
	Hydranal 10.0 (B18F011-SRM2)_D0505ACRLKTF_20180612-122659.mdet	6/12/2018 12:27 PM	MDET File	78 KB		
	Hydranal 10.0 (B18F011-SRM3)_D0505ACRLKTF_20180612-122659.mdet	6/12/2018 12:27 PM	MDET File	78 KB		
	Hydranal 10.0 (B18F011-SRM4)_D0505ACRLKTF_20180612-122659.mdet	6/12/2018 12:27 PM	MDET File	78 KB		
	Water_D0505ACRLKTF_20180612-122700.mdet	6/12/2018 12:27 PM	MDET File	227 KB		

DD)



## UNITED STATES ENVIRONMENTAL PROTECTION AGENCY

## REGION 5 CHICAGO REGIONAL LABORATORY

536 SOUTH CLARK STREET

CHICAGO, ILLINOIS 60605

**Data Verification Checklist****Work Order:** 1806004**Analysis:** Ignitability by Setaflash**Project:** Milliken Chemical Dewey Plant**CRL SOP Used:** AIG048A Version 5**Electronic Pathway(s):**

I:\Projects (New File Structure)\Milliken Chemical Dewey Plant\1. A&amp;I\CBreslin\1806004\Setaflash 3\flashpoint

**File Name(s):**

Not Applicable

**General Information**

Question #	Question	Analyst Response (YES/ NO/ NA)	Reviewer Response (Circle Response)
1	Did samples meet the laboratory's standard conditions of sample acceptability upon receipt?	YES	<input checked="" type="checkbox"/> YES NO / NA
2	Was customer contact communication included?	YES	<input checked="" type="checkbox"/> YES NO / NA
3	Are all field sample ID numbers cross-referenced to the laboratory ID numbers?	YES	<input checked="" type="checkbox"/> YES NO / NA
4	Were all samples prepared and analyzed within holding times?	NA	YES / NO / <input checked="" type="checkbox"/> NA

**Sample Results**

Question #	Question	Analyst Response (YES/ NO/ NA)	Reviewer Response (Circle Response)
5	Were calculations checked?	YES	<input checked="" type="checkbox"/> YES NO / NA
6	Were all analyte results checked?	YES	<input checked="" type="checkbox"/> YES NO / NA
7	Were all results for soil and sediment samples reported on a dry weight basis?	NA	YES / NO / <input checked="" type="checkbox"/> NA
8	Were % moisture (or solids) reported for all soil and sediment samples?	NA	YES / NO / <input checked="" type="checkbox"/> NA
9	Other than those results < RL, were all other raw values bracketed by calibration standards?	YES	<input checked="" type="checkbox"/> YES NO / NA
10	Are the RLs for each method analyte included in the laboratory data package?	NA	YES / NO / <input checked="" type="checkbox"/> NA
11	Are MDLs/RLs adjusted for dilutions?	NA	YES / NO / <input checked="" type="checkbox"/> NA
12	Were the raw data (for example, chromatograms, spectral data) reviewed?	YES	<input checked="" type="checkbox"/> YES NO / NA
13	Were data associated with manual integrations flagged on the raw data?	NA	YES / NO / <input checked="" type="checkbox"/> NA

**Standards**

Question #	Question	Analyst Response (YES/ NO/ NA)	Reviewer Response (Circle Response)
14	Are all standards used in the analyses NIST-traceable or obtained from other appropriate sources as specified in the analytical SOP?	YES	<input checked="" type="checkbox"/> YES NO / NA
15	Were standard/reagent preparations checked (preparation date, expiration date, parent standard IDs, etc.), if applicable?	YES	<input checked="" type="checkbox"/> YES NO / NA

**Batch Quality Control**

Question #	Question	Analyst Response (YES/ NO/ NA)	Reviewer Response (Circle Response)
16	Were appropriate type(s) of blanks analyzed?	NA	YES / NO / <input checked="" type="checkbox"/> NA



### Batch Quality Control

Question #	Question	Analyst Response (YES/ NO/ NA)	Reviewer Response (Circle Response)
17	Were blanks analyzed at the appropriate frequency?	NA	YES / NO <input checked="" type="radio"/> NA
18	Were method blanks taken through the entire analytical process, including preparation and, if applicable, cleanup procedures?	NA	YES / NO <input checked="" type="radio"/> NA
19	Were blank concentrations <= MDL or RL, as applicable in the analytical SOP?	NA	YES / NO <input checked="" type="radio"/> NA
20	Was each LCS/SRM taken through the entire analytical procedure, including prep and cleanup steps?	YES	<input checked="" type="radio"/> YES / NO / NA
21	Were LCSs/SRMs analyzed at the required frequency?	YES	<input checked="" type="radio"/> YES / NO / NA
22	Were LCS/SRM %Rs within the laboratory QC limits or other acceptance criteria?	YES	<input checked="" type="radio"/> YES / NO / NA
23	Were the project/method specified analytes included in the MS?	NA	YES / NO <input checked="" type="radio"/> NA
24	Were MS analyzed at the appropriate frequency?	NA	YES / NO <input checked="" type="radio"/> NA
25	Were MS %Rs within the laboratory QC limits?	NA	YES / NO <input checked="" type="radio"/> NA
26	Were appropriate analytical duplicates analyzed for each matrix?	YES	<input checked="" type="radio"/> YES / NO / NA
27	Were analytical duplicates analyzed at the appropriate frequency?	YES	<input checked="" type="radio"/> YES / NO / NA
28	Were RPDs or relative standard deviations within the laboratory QC limits or other acceptance criteria?	YES	<input checked="" type="radio"/> YES / NO / NA
29	Were RLs analyzed at the appropriate frequency?	NA	YES / NO <input checked="" type="radio"/> NA
30	Were RL recoveries within the laboratory QC limits?	NA	YES / NO <input checked="" type="radio"/> NA

### Calibration

Question #	Question	Analyst Response (YES/ NO/ NA)	Reviewer Response (Circle Response)
31	Were calibration correlation coefficient or other acceptance criteria met?	YES	<input checked="" type="radio"/> YES / NO / NA
32	Was the number of calibration standards recommended in the method used for all analytes?	NA	YES / NO <input checked="" type="radio"/> NA
33	Were all points generated between the lowest and highest standard used to calculate the curve?	NA	YES / NO <input checked="" type="radio"/> NA
34	Are calibration data available for all instruments used?	YES	<input checked="" type="radio"/> YES / NO / NA

### Calibration Verification

Question #	Question	Analyst Response (YES/ NO/ NA)	Reviewer Response (Circle Response)
35	Was the absolute value of the analyte concentration in the ICB/CCB <= MDL or RL, as applicable in the analytical SOP?	NA	YES / NO <input checked="" type="radio"/> NA
36	Was the calibration curve verified for each analyte?	NA	YES / NO <input checked="" type="radio"/> NA
37	Has the calibration curve been verified using an appropriate second source standard?	NA	YES / NO <input checked="" type="radio"/> NA
38	Were ICV/CCV analyzed at the method-required frequency?	NA	YES / NO <input checked="" type="radio"/> NA
39	Were ICV/CCV %R within the laboratory QC limits?	NA	YES / NO <input checked="" type="radio"/> NA

### Quality Control

Question #	Question	Analyst Response (YES/ NO/ NA)	Reviewer Response (Circle Response)
40	QC limits in LIMS checked against the SOP?	YES	<input checked="" type="radio"/> YES / NO / NA
41	Were percent differences, recoveries, and the linearity within the QC limits specified in the method?	YES	<input checked="" type="radio"/> YES / NO / NA

### Quality Control

Question #	Question	Analyst Response (YES/ NO/ NA)	Reviewer Response (Circle Response)
42	Were QC charts checked if a QC audit, excluding matrix QC (MS or DUP), was out of limit?	NA	YES / NO / <b>NA</b>

### Supporting Data

Question #	Question	Analyst Response (YES/ NO/ NA)	Reviewer Response (Circle Response)
43	Were support equipment data (balance verification, data logs, logbook entries, etc.) included in data package?	YES	<b>YES</b> / NO / NA
44	Were approved spreadsheet(s) used?	YES	<b>YES</b> / NO / NA

### Document Verification

Question #	Question	Analyst Response (YES/ NO/ NA)	Reviewer Response (Circle Response)
45	Is the MDL or RL study up-to-date for each reported analyte?	NA	YES / NO / <b>NA</b>
46	Is documentation of the analyst's capability up-to-date and on file? <b>ID 14749</b>	YES	<b>YES</b> / NO / NA
47	Are the procedures for compound/analyte identification documented?	YES	<b>YES</b> / NO / NA
48	Are laboratory SOPs current and on file for the method performed?	YES	<b>YES</b> / NO / NA
49	Were data evaluated against project QAPP or Sample Plan and documented in case narrative?	YES	<b>YES</b> / NO / NA

### Accreditation

Question #	Question	Analyst Response (YES/ NO/ NA)	Reviewer Response (Circle Response)
50	ANAB logo used appropriately?	YES	<b>YES</b> / NO / NA

Signature and Date:

**REVIEWED**

By Colin Breslin at 9:33 am, Jul 17, 2018

**REVIEWED**

By Francis A Awanya at 1:05 pm, Jul 18, 2018



## UNITED STATES ENVIRONMENTAL PROTECTION AGENCY

## REGION 5 CHICAGO REGIONAL LABORATORY

536 SOUTH CLARK STREET

CHICAGO, ILLINOIS 60605

**Data Verification Checklist****Work Order:** 1806004**Analysis:** Water content, Karl Fisher Titration**Project:** Milliken Chemical Dewey Plant**CRL SOP Used:** AIG015A Version 4**Electronic Pathway(s):**

I:\Projects (New File Structure)\Milliken Chemical Dewey Plant\1. A&amp;I\CBreslin\1806004\KFT\Water Content

**File Name(s):**

See Memo-to-File: Data Archive Paths

**General Information**

Question #	Question	Analyst Response (YES/ NO/ NA)	Reviewer Response (Circle Response)
1	Did samples meet the laboratory's standard conditions of sample acceptability upon receipt?	YES	<input checked="" type="checkbox"/> YES / NO / NA
2	Was customer contact communication included?	YES	<input checked="" type="checkbox"/> YES / NO / NA
3	Are all field sample ID numbers cross-referenced to the laboratory ID numbers?	YES	<input checked="" type="checkbox"/> YES / NO / NA
4	Were all samples prepared and analyzed within holding times?	NA	YES / NO / <input checked="" type="checkbox"/> NA

**Sample Results**

Question #	Question	Analyst Response (YES/ NO/ NA)	Reviewer Response (Circle Response)
5	Were calculations checked?	YES	<input checked="" type="checkbox"/> YES / NO / NA
6	Were all analyte results checked?	YES	<input checked="" type="checkbox"/> YES / NO / NA
7	Were all results for soil and sediment samples reported on a dry weight basis?	NA	YES / NO / <input checked="" type="checkbox"/> NA
8	Were % moisture (or solids) reported for all soil and sediment samples?	NA	YES / NO / <input checked="" type="checkbox"/> NA
9	Other than those results < RL, were all other raw values bracketed by calibration standards?	YES	<input checked="" type="checkbox"/> YES / NO / NA
10	Are the RLs for each method analyte included in the laboratory data package?	YES	<input checked="" type="checkbox"/> YES / NO / NA
11	Are MDLs/RLs adjusted for dilutions?	NA	YES / NO / <input checked="" type="checkbox"/> NA
12	Were the raw data (for example, chromatograms, spectral data) reviewed?	YES	<input checked="" type="checkbox"/> YES / NO / NA
13	Were data associated with manual integrations flagged on the raw data?	NA	YES / NO / <input checked="" type="checkbox"/> NA

**Standards**

Question #	Question	Analyst Response (YES/ NO/ NA)	Reviewer Response (Circle Response)
14	Are all standards used in the analyses NIST-traceable or obtained from other appropriate sources as specified in the analytical SOP?	YES	<input checked="" type="checkbox"/> YES / NO / NA
15	Were standard/reagent preparations checked (preparation date, expiration date, parent standard IDs, etc.), if applicable?	YES	<input checked="" type="checkbox"/> YES / NO / NA

**Batch Quality Control**

Question #	Question	Analyst Response (YES/ NO/ NA)	Reviewer Response (Circle Response)
16	Were appropriate type(s) of blanks analyzed?	NA	YES / NO / <input checked="" type="checkbox"/> NA

### Batch Quality Control

Question #	Question	Analyst Response (YES/ NO/ NA)	Reviewer Response (Circle Response)
17	Were blanks analyzed at the appropriate frequency?	NA	YES / NO / <input checked="" type="radio"/> NA
18	Were method blanks taken through the entire analytical process, including preparation and, if applicable, cleanup procedures?	NA	YES / NO / <input checked="" type="radio"/> NA
19	Were blank concentrations <= MDL or RL, as applicable in the analytical SOP?	NA	YES / NO / <input checked="" type="radio"/> NA
20	Was each LCS/SRM taken through the entire analytical procedure, including prep and cleanup steps?	YES	<input checked="" type="radio"/> YES / NO / NA
21	Were LCSs/SRMs analyzed at the required frequency?	YES	<input checked="" type="radio"/> YES / NO / NA
22	Were LCS/SRM %Rs within the laboratory QC limits or other acceptance criteria?	YES	<input checked="" type="radio"/> YES / NO / NA
23	Were the project/method specified analytes included in the MS?	NA	YES / NO / <input checked="" type="radio"/> NA
24	Were MS analyzed at the appropriate frequency?	NA	YES / NO / <input checked="" type="radio"/> NA
25	Were MS %Rs within the laboratory QC limits?	NA	YES / NO / <input checked="" type="radio"/> NA
26	Were appropriate analytical duplicates analyzed for each matrix?	YES	<input checked="" type="radio"/> YES / NO / NA
27	Were analytical duplicates analyzed at the appropriate frequency?	YES	<input checked="" type="radio"/> YES / NO / NA
28	Were RPDs or relative standard deviations within the laboratory QC limits or other acceptance criteria?	YES	<input checked="" type="radio"/> YES / NO / NA
29	Were RLs analyzed at the appropriate frequency?	NA	YES / NO / <input checked="" type="radio"/> NA
30	Were RL recoveries within the laboratory QC limits?	NA	YES / NO / <input checked="" type="radio"/> NA

### Calibration

Question #	Question	Analyst Response (YES/ NO/ NA)	Reviewer Response (Circle Response)
31	Were calibration correlation coefficient or other acceptance criteria met?	YES	<input checked="" type="radio"/> YES / NO / NA
32	Was the number of calibration standards recommended in the method used for all analytes?	NA	YES / NO / <input checked="" type="radio"/> NA
33	Were all points generated between the lowest and highest standard used to calculate the curve?	NA	YES / NO / <input checked="" type="radio"/> NA
34	Are calibration data available for all instruments used?	YES	<input checked="" type="radio"/> YES / NO / NA

### Calibration Verification

Question #	Question	Analyst Response (YES/ NO/ NA)	Reviewer Response (Circle Response)
35	Was the absolute value of the analyte concentration in the ICB/CCB <= MDL or RL, as applicable in the analytical SOP?	NA	YES / NO / <input checked="" type="radio"/> NA
36	Was the calibration curve verified for each analyte?	NA	YES / NO / <input checked="" type="radio"/> NA
37	Has the calibration curve been verified using an appropriate second source standard?	NA	YES / NO / <input checked="" type="radio"/> NA
38	Were ICV/CCV analyzed at the method-required frequency?	NA	YES / NO / <input checked="" type="radio"/> NA
39	Were ICV/CCV %R within the laboratory QC limits?	NA	YES / NO / <input checked="" type="radio"/> NA

### Quality Control

Question #	Question	Analyst Response (YES/ NO/ NA)	Reviewer Response (Circle Response)
40	QC limits in LIMS checked against the SOP?	YES	<input checked="" type="radio"/> YES / NO / NA
41	Were percent differences, recoveries, and the linearity within the QC limits specified in the method?	YES	<input checked="" type="radio"/> YES / NO / NA

### Quality Control

Question #	Question	Analyst Response (YES/ NO/ NA)	Reviewer Response (Circle Response)
42	Were QC charts checked if a QC audit, excluding matrix QC (MS or DUP), was out of limit?	NA	YES / NO / <b>NA</b>

### Supporting Data

Question #	Question	Analyst Response (YES/ NO/ NA)	Reviewer Response (Circle Response)
43	Were support equipment data (balance verification, data logs, logbook entries, etc.) included in data package?	YES	<b>YES</b> / NO / NA
44	Were approved spreadsheet(s) used?	YES	<b>YES</b> / NO / NA

### Document Verification

Question #	Question	Analyst Response (YES/ NO/ NA)	Reviewer Response (Circle Response)
45	Is the MDL or RL study up-to-date for each reported analyte?	NA	YES / NO / <b>NA</b>
46	Is documentation of the analyst's capability up-to-date and on file? <b>ID 15096</b>	YES	<b>YES</b> / NO / NA
47	Are the procedures for compound/analyte identification documented?	YES	<b>YES</b> / NO / NA
48	Are laboratory SOPs current and on file for the method performed?	YES	<b>YES</b> / NO / NA
49	Were data evaluated against project QAPP or Sample Plan and documented in case narrative?	YES	<b>YES</b> / NO / NA

### Accreditation

Question #	Question	Analyst Response (YES/ NO/ NA)	Reviewer Response (Circle Response)
50	ANAB logo used appropriately?	YES	<b>YES</b> / NO / NA

Signature and Date:

**REVIEWED**  
By Colin Breslin at 10:38 am, Jul 17, 2018

**REVIEWED**  
By Francis A Awanya at 1:06 pm, Jul 18, 2018



Environmental Protection Agency Region 5  
**Chicago Regional Laboratory**

Date: 6/13/2018

536 South Clark Street, Chicago, IL 60605  
Phone:(312)353-8370 Fax:(312)886-2591

**AUTOMATED ANALYZER BENCH SHEET**

**Batch Number: B18F004**

**Analyses Included On This Benchsheet**

**Ignitability by Setaflash**

Date: 6/13/2018

Analyst (initials): CB

Sample Number	Source ID	Spike Type	Spike ID	Spike Amount (uL)	Final Volume (mL)	Dilution (mL)		Comments
						Initial	Final	
1805008-01								pH = 4
1805008-02								pH = 4
1805008-03								pH = 7
1805008-04								pH = 7
1805008-05								pH = 4
1805008-06								pH = 9
1806004-01								pH = 4
1806004-02								pH = 4
B18F004-DUP1	1805008-02							
B18F004-DUP2	1806004-01							
B18F004-SRM1		Static	17E2502		2			
B18F004-SRM2		Static	17E2502		2			
B18F004-SRM3		Static	17E2502		2			avg of SRM1 & SRM2

Batch Comments:



Environmental Protection Agency Region 5  
**Chicago Regional Laboratory**

Date: 6/13/2018

536 South Clark Street, Chicago, IL 60605  
Phone:(312)353-8370 Fax:(312)886-2591

**AUTOMATED ANALYZER BENCH SHEET**

**Batch Number: B18F004**

**Analyses Included On This Benchsheet**

**Ignitability by Setaflash**

Date: 6/13/2018

Analyst (initials): CB

Sample Number	Source ID	Spike Type	Spike ID	Spike Amount (uL)	Final Volume (mL)	Dilution (mL)		Comments
						Initial	Final	
B18F004-SRM4		Static	17E2502		2			
B18F004-SRM5		Static	17E2502		2			
B18F004-SRM6		Static	17E2502		2			avg of SRM4 & SRM5
B18F004-SRM7		Static	17E2502		2			
B18F004-SRM8		Static	17E2502		2			
B18F004-SRM9		Static	17E2502		2			avg of SRM7 & SRM8

**PREPARATION REAGENTS/STANDARDS/PIPETTES:**

**17E2502:** PARAXYLENE FLASH POINT CHECK

**Prepared: May-25-17**

**Expires: May-25-27**

**17C1307:** Barometer - Digital

**Prepared: Mar-13-17**

**Expires: Feb-03-19**

Batch Comments:

## Working Bench Sheet

Work Order(s): 1805008  
LIMS Batch ID: B18F004

Date: 6/6/18

Analyst: CB

Sample ID: SRM, LIMS # 17E2502  
B18F004-SRM1

Sample ID: SRM, LIMS # 17E2502  
B18F004 - SRM2

[illegible][illegible]

Sample ID: 1805008-C1

### **Preliminary (Range Finding) -**

[illegible]

Sample ID: 1805008-01

**Finite Flash Point -**

[illegible]

no finite.  
No Flash  
observed  
during  
range  
finding.  
-CB6/6/18



## Flash Point by AIG048A

Continuation from previous page

Date: 6/6/18Analyst: CBSample ID: 1501805008-02**Preliminary (Range Finding) -**

Temperature (°F)	Flash (Y/N)	Pressure (specify units)
78	NO	990 mb
87	NO	990 mb
96	NO	990 mb
105	NO	990 mb
114	NO	990 mb
123	NO	991 mb
132	NO	991 mb
141	NO	991 mb
149	YES	991 mb

Sample ID: \_\_\_\_\_

**Preliminary (Range Finding) -**

Temperature (°F)	Flash (Y/N)	Pressure (specify units)

Range Finding NA

For duplicate.

-CB 6/6/18

Sample ID: 1805008-02**Finite Flash Point -**

Temperature (°F)	Flash (Y/N)	Pressure (specify units)
141	NO	991 mb
142	NO	991 mb
143	NO	991 mb
144	NO	991 mb
145	YES	991 mb
143	NO	991 mb
144	NO	990 mb
145	NO	990 mb
146	YES	990 mb

Sample ID: B18F004-D9P1 (Source = 1805008-02)**Finite Flash Point -**

Temperature (°F)	Flash (Y/N)	Pressure (specify units)
141	NO	989 mb
142	NO	989 mb
143	NO	989 mb
144	NO	988 mb
145	NO	988 mb
146	NO	988 mb
147	YES	988 mb
144	NO	988 mb
145	NO	989 mb
146	NO	989 mb
147	YES	988 mb

## Working Bench Sheet

Analyst: CB

Sample ID: SRM, LIMS # 17E2502  
B18F004-SRM5

[illegible]

86/7/18

Sample ID: 1805008-03

No Finite  
determination.  
-CB 6/7/18

[illegible]

red lower layer. black top layer too thin to phase separate. Able to mix and emulsify top layer in bottom layer to take a representative sample aliquot.

Flash Point by AIG048A

Continuation from previous page

Date: 6/7/18

Analyst: CB

Sample ID: 1805008-04

Preliminary (Range Finding) -

CB 6/7/18

Temperature (°F)	Flash (Y/N)	Pressure (specify units)
79.78	NO	992 mb
87	NO	992 mb
96	NO	992 mb
105	NO	992 mb
114	NO	992 mb
123	NO	992 mb
132	NO	992 mb
141	NO	992 mb
149	NO	992 mb
140	NO	992 mb
149	NO	992 mb

Sample ID: 1805008-04

Finite Flash Point -

Temperature (°F)	Flash (Y/N)	Pressure (specify units)

no finite.  
No Flash during range Finding.  
-CB 6/7/18

Sample ID: 1805008-05

Preliminary (Range Finding) -

Temperature (°F)	Flash (Y/N)	Pressure (specify units)
81	NO	993 mb
90	NO	993 mb
99	NO	993 mb
108	NO	993 mb
117	NO	992 mb
126	NO	992 mb
135	NO	992 mb
144	NO	992 mb
149	NO	992 mb
140	NO	992 mb
149	NO	992 mb

Sample ID: 1805008-05

Finite Flash Point -

Temperature (°F)	Flash (Y/N)	Pressure (specify units)

No finite.  
No Flash during range Finding.  
-CB 6/7/18



## Working Bench Sheet

Analyst: CB

Sample ID: SRM, LIMS # 17E2502  
BIGFOOT-SRM8

[illegible][illegible]

Sample ID: 1405008-06

**Finite Flash Point -**

[illegible][illegible]

no finite because  
no flash during range  
Finding. - CB 6/8/18

## Continuation from previous page

Analyst: CB

### **Preliminary (Range Finding) -**

[illegible]

### Preliminary (Range Finding) -

[illegible]

Range Finding NA for  
duplicate. Finite determination  
duplicated. -CB 6/8/18

**Finite Flash Point -**

[illegible]

996 mb CB 6/8/18

**Finite Flash Point -**

[illegible]

1806004-0

## Continuation from previous page

Analyst: CB

### **Preliminary (Range Finding) -**

[illegible]

### Preliminary (Range Finding) -

[illegible]

**Finite Flash Point -**

[illegible]

**Finite Flash Point -**

[illegible]



Environmental Protection Agency Region 5  
**Chicago Regional Laboratory**

Date: 6/13/2018

536 South Clark Street, Chicago, IL 60605  
Phone:(312)353-8370 Fax:(312)886-2591

**AUTOMATED ANALYZER BENCH SHEET**

**Batch Number: B18F011**

**Analyses Included On This Benchsheet**

**Water content,Karl Fisher Titration**

Date: 6/13/2018

Analyst (initials): CB

Sample Number	Source ID	Spike Type	Spike ID	Spike Amount (uL)	Final Volume (mL)	Dilution (mL)		Comments
						Initial	Final	
1805008-01								pH = 4
1805008-02								pH = 4
1805008-03								pH = 7
1805008-04								pH = 7
1805008-05								pH = 4
1805008-06								pH = 9
1806004-01								pH = 4
1806004-02								pH = 4
1806009-01								same as B18F011-SRM1
1806009-02								same as B18F011-SRM2
1806009-03								same as B18F011-SRM3
1806009-04								same as B18F011-SRM4
B18F011-DUP1	1805008-01							

Batch Comments: updated with P&A WO 1806009.



Environmental Protection Agency Region 5  
**Chicago Regional Laboratory**

Date: 6/13/2018

536 South Clark Street, Chicago, IL 60605  
Phone:(312)353-8370 Fax:(312)886-2591

**AUTOMATED ANALYZER BENCH SHEET**

**Batch Number: B18F011**

**Analyses Included On This Benchsheet**

**Water content,Karl Fisher Titration**

Date: 6/13/2018

Analyst (initials): CB

Sample Number	Source ID	Spike Type	Spike ID	Spike Amount (uL)	Final Volume (mL)	Dilution (mL)		Comments
						Initial	Final	
B18F011-DUP2	1805008-02							
B18F011-DUP3	1805008-03							
B18F011-DUP4	1805008-04							
B18F011-DUP5	1805008-05							
B18F011-DUP6	1805008-06							
B18F011-DUP7	1806004-01							
B18F011-DUP8	1806004-02							
B18F011-SRM1		Static	17E1709		0.05			P&A1
B18F011-SRM2		Static	17E1709		0.05			P&A2
B18F011-SRM3		Static	17E1709		0.05			P&A3
B18F011-SRM4		Static	17E1709		0.05			P&A4

Batch Comments: updated with P&A WO 1806009.





Environmental Protection Agency Region 5  
**Chicago Regional Laboratory**

Date: 6/13/2018

536 South Clark Street, Chicago, IL 60605  
Phone:(312)353-8370 Fax:(312)886-2591

**AUTOMATED ANALYZER BENCH SHEET**

**Batch Number: B18F011**

**Analyses Included On This Benchsheet**

**Water content,Karl Fisher Titration**

Date: 6/13/2018

Analyst (initials): CB

Sample Number	Source ID	Spike Type	Spike ID	Spike Amount (uL)	Final Volume (mL)	Dilution (mL)		Comments
						Initial	Final	
<b><u>PREPARATION REAGENTS/STANDARDS/PIPETTES:</u></b>								
<b><u>17E1709:</u></b>	Water Content - Hydranal Water Standard 10.0			Prepared: Jan-19-16		Expires: Dec-23-20		
<b><u>17E1710:</u></b>	HYDRANAL - Composite 5 K			Prepared: Jul-07-16		Expires: Jun-22-19		
<b><u>18D1804:</u></b>	HYDRANAL - Working Medium K			Prepared: Oct-04-17		Expires: Sep-24-19		

Batch Comments: updated with P&A WO 1806009.



Environmental Protection Agency Region 5  
**Chicago Regional Laboratory**

536 South Clark Street, Chicago, IL 60605  
Phone: (312) 353-8370 Fax: (312) 886-2591

Date: 6/12/2018

Water Content by Karl Fischer Titration

Date: 6/12/18

Analyst (initials): CB

**Batch Number: B18F011**

**Analyses Included On This Benchsheet**

**Water content, Karl Fisher Titration**

Sample number (LIMS ID)	Duplicate Source sample (LIMS ID)	Initial syringe weight (g)	Final syringe weight (g)	Sample weight (g)	Comments
1805008-01				0.05	pH = 4
1805008-02				0.05	pH = 4
1805008-03				0.05	pH = 7
1805008-04				0.05	pH = 7 able to emulsify top layer for representative sample
1805008-05				0.05	pH = 4
1805008-06				0.05	pH = 9
1806004-01				0.05	pH = 4
1806004-02				0.05	pH = 4
B18F011-DUP1	1805008-01			0.05	

Batch Comments:

Comments:

recorded in instrument  
software.  
-CB 6/12/18



Environmental Protection Agency Region 5  
**Chicago Regional Laboratory**

Date: 6/12/2018

536 South Clark Street, Chicago, IL 60605  
Phone: (312) 353-8370 Fax: (312) 886-2591

**Water Content by Karl Fischer Titration**

Date: 6/12/18  
Analyst (initials): CB

**Batch Number: B18F011**  
**Analyses Included On This Benchsheet**  
**Water content, Karl Fisher Titration**

Sample number (LIMS ID)	Duplicate Source sample (LIMS ID)	Initial syringe weight (g)	Final syringe weight (g)	Sample weight (g)	Comments
B18F011-DUP2	1805008-02			0.05	
B18F011-DUP3	1805008-03			0.05	
B18F011-DUP4	1805008-04			0.05	
B18F011-DUP5	1805008-05			0.05	
B18F011-DUP6	1805008-06			0.05	
B18F011-DUP7	1806004-01			0.05	
B18F011-DUP8	1806004-02			0.05	
B18F011-SRM1	-			0.05	P&A1
B18F011-SRM2	-			0.05	P&A2

Batch Comments:

Comments:

CB 6/12/18  
recorded in software.



Environmental Protection Agency Region 5  
**Chicago Regional Laboratory**

Date: 6/12/2018

536 South Clark Street, Chicago, IL 60605  
Phone: (312) 353-8370 Fax: (312) 886-2591

**Water Content by Karl Fischer Titration**

Date: 6/12/18  
Analyst (initials): CB

**Batch Number: B18F011**  
**Analyses Included On This Benchsheet**  
**Water content, Karl Fisher Titration**

Sample number (LIMS ID)	Duplicate Source sample (LIMS ID)	Initial syringe weight (g)	Final syringe weight (g)	Sample weight (g)	Comments
B18F011-SRM3	-			<del>0.05</del>	P&A3
B18F011-SRM4	-			<del>0.05</del>	P&A4

**PREPARATION REAGENTS/STANDARDS/PIPETTES:**

**17E1709:** Water Content - Hydranal Water  
Standard 10.0  
**17E1710:** HYDRANAL - Composite 5 K  
**18D1804:** HYDRANAL - Working Medium K

CB 6/12/18  
recorded in instrument  
software.

Batch Comments:

Comments:



	<u>Standard Weight</u>	<u>Actual Weight</u>
12/28/17	0.1g	0.1000g
CB	0.5g	0.5001g
Weight Set	1g	1.0000g
Troemner	5g	4.9999g
S/N 4000017548	10g	10.0001g
	50g	49.9999g
	100g	100.0000g

CB 1/9/18 • Balance OK after building water damage.

	<u>Standard Weight</u>	<u>Reading</u>
CB 5/17/18	0.1g	0.0999g
Weight Set	0.5g	0.5000g
Troemner	1g	1.0000g
S/N 4000020557	5g	5.0000g
	10g	9.9999g
	50g	50.0001g
	100g	100.0001g

CB 6/6/18	0.1g	0.1000g
Weight Set	0.5g	0.5000g
Troemner	1g	1.0001g
S/N 4000020557	5g	5.0000g
	10g	10.0000g
	50g	49.9999g
	100g	99.9998g

22

	<u>Standard Weight</u>	<u>Actual Weight</u>
CB 6/7/18	0.1g	0.1000g
Weight Set	0.5g	0.4999g
Troemner	1g	0.9999g
S/N 4000020557	5g	5.0000g
	10g	10.0000g
	50g	50.0000g
	100g	100.0000g

CB 6/8/18	0.1g	0.1000g
CB 6/8/18	0.5g	0.5000g
Weight Set	1g	0.9999g
Troemner	5g	5.0001g
S/N 40000 20557	10g	10.0001g
	50g	50.0000g
	100g	99.9999g

CB 6/12/18	0.1g	0.1000g
Weight set	0.5g	0.4999g
Troemner	1g	1.0001g
S/N 40000 20557	5g	5.0000g
	10g	10.0000g
	50g	50.0001g
	100g	99.9999g



- CB  
6/14/17 • Analyzed B17F023.  
• System ran well, only issues with difficult sample matrices were encountered. Titer strength, conditioning response, and Hydram 10.0 (LIMS 17E1709) all performed well and/or met QC limits.
- CB  
6/15/17 • Analyzed B17F024  
• Titer strength determination took multiple attempts. Cleaned apparatus and repeated until passing.
- CB  
7/26/17 • Analyzed B17G047  
• System ran well. Conditioning was good.
- CB  
12/28/17 • Pump to Fill titration vessel was not working.  
Opened 803 Ti Stand housing and Tightened down switch. Filling pump began working.  
• Analyzed B17L033 CB 12/28/17  
• Titer strength determination required two attempts. Afterwards conditioning was successful and run was OK.
- CB  
1/9/18 • Instrument and computer OK after building water damage. Balance #21 OK too.
- CB  
6/12/18 • Analyzed B18F011. System Ran Well.

6

Setaflash 3 instrument logbook

5/17/18 • Re-installed Thermometer (S/N 85C020 exp 9/15/18)  
CB after being out for calibration.

• Syringe check. Syringe ID = N-1244.

Balance #21. 1.9414 g. OK at  $2.0 \pm 0.1$  mL.

• Barometer ID = LIMS 17C1307.

6/6/18 • Syringe ID = N-1244

CB

• Syringe check on Balance #21, 1.9518 g.  
OK at  $2.0 \pm 0.1$  mL.

• Barometer ID = LIMS 17C1307

6/7/18

CB

• Syringe ID = N-1244

• Balance #21

• 1.9272 g  $\rightarrow$  OK at  $2.0 \pm 0.1$  mL.

• Barometer ID = LIMS 17C1307

6/8/18

CB

• Syringe ID = N-1244

• Balance #21

• 1.9420 g  $\rightarrow$  OK at  $2.0 \pm 0.1$  mL.

• Barometer ID = LIMS 17C1307



## ANALYSIS SEQUENCE

Printed: 6/13/2018 11:40:30AM

18F1302

Instrument: KFTAnalysis Date: Jun-12-18Calibration ID: UNASSIGNEDMatrix: Water

Order	Container	Lab Number	Sample name or QC sample designation	Analysis	STD ID	Comments
1		B18F011-SRM1	Reference	QC		
2	A	1806009-01	P&A1	Water content,Karl Fisher Titration		
3		B18F011-SRM2	Reference	QC		
4	A	1806009-02	P&A2	Water content,Karl Fisher Titration		
5		B18F011-SRM3	Reference	QC		
6	A	1806009-03	P&A3	Water content,Karl Fisher Titration		
7		B18F011-SRM4	Reference	QC		
8	A	1806009-04	P&A4	Water content,Karl Fisher Titration		
9	C	1805008-02	HCC02	Water content,Karl Fisher Titration		pH = 4
10		B18F011-DUP2	Duplicate	QC		
11	B	1806004-01	01 Tetramer	Water content,Karl Fisher Titration		pH = 4
12		B18F011-DUP7	Duplicate	QC		
13	B	1806004-02	02 Tetramer	Water content,Karl Fisher Titration		pH = 4
14		B18F011-DUP8	Duplicate	QC		
15	B	1805008-01	HCC01	Water content,Karl Fisher Titration		pH = 4
16		B18F011-DUP1	Duplicate	QC		
17	B	1805008-03	HCC03	Water content,Karl Fisher Titration		pH = 7
18		B18F011-DUP3	Duplicate	QC		
19	C	1805008-04	HCC04	Water content,Karl Fisher Titration		pH = 7
20		B18F011-DUP4	Duplicate	QC		
21	C	1805008-05	HCC05	Water content,Karl Fisher Titration		pH = 4
22		B18F011-DUP5	Duplicate	QC		

Order	Container	Lab Number	Sample name <u>or</u> QC sample designation	Analysis	STD ID	Comments
23	B	1805008-06	HCC06	Water content,Karl Fisher Titration		pH = 9
24		B18F011-DUP6	Duplicate	QC		

**PREPARATION REAGENTS/STANDARDS/PIPETTES:**

**17E1709:**    Water Content - Hydranal Water Standard 10.0

**Prepared: Jan-19-16**

**Expires: Dec-23-20**

***REVIEWED***

***By Colin Breslin at 11:41 am, Jun 13, 2018***

---

Analyst Initials and Date

Determination overview

	Determination start	Sample ID	Method name	Titer (mg/mL)	KF Titer Mean value	Titer s(rei)	% Water	mg Water	Sample size (g)	mL to EP	Remarks
1	2018-06-12 12:20:35 UTC-5	B18F011-DUP6	KF SAMPLE				96.29	962.94	-0.0512	10.25	
2	2018-06-12 12:17:51 UTC-5	1805008-06	KF SAMPLE				96.06	960.57	-0.0606	12.10	
3	2018-06-12 12:13:33 UTC-5	B18F011-DUP5	KF SAMPLE				0.05	0.50	-4.2541	0.44	
4	2018-06-12 12:10:12 UTC-5	1805008-05	KF SAMPLE				0.05	0.50	-4.0966	0.42	
5	2018-06-12 12:06:36 UTC-5	1806008-06	KF SAMPLE				0.26	2.67	-0.0412	0.09	
6	2018-06-12 12:01:26 UTC-5	B18F011-DUP4	KF SAMPLE				55.34	553.43	-0.0591	6.80	1
7	2018-06-12 11:56:50 UTC-5	1805008-04	KF SAMPLE				55.58	555.77	-0.0615	7.11	
8	2018-06-12 11:52:39 UTC-5	B18F011-DUP3	KF SAMPLE				96.84	968.37	-0.0423	8.52	
9	2018-06-12 11:49:34 UTC-5	1805008-03	KF SAMPLE				96.49	964.94	-0.0425	8.53	
10	2018-06-12 11:43:39 UTC-5	1806008-08	KF SAMPLE				invalid	invalid	-1.0466	invalid	
11	2018-06-12 11:40:34 UTC-5	B18F011-DUP1	KF SAMPLE				0.01	0.12	-4.2388	0.10	2
12	2018-06-12 11:37:21 UTC-5	1805008-01	KF SAMPLE				0.01	0.10	-3.9470	0.08	
13	2018-06-12 11:34:18 UTC-5	1806008-04	KF SAMPLE				0.02	0.18	-0.8232	0.09	
14	2018-06-12 11:26:15 UTC-5	B18F011-DUP8	KF SAMPLE				0.01	0.09	-3.9783	0.07	3
15	2018-06-12 11:23:07 UTC-5	1806004-02	KF SAMPLE				0.01	0.08	-3.8427	0.06	
16	2018-06-12 11:18:23 UTC-5	B18F011-DUP7	KF SAMPLE				0.01	0.10	-3.8794	0.08	
17	2018-06-12 11:14:52 UTC-5	1806004-01	KF SAMPLE				0.01	0.09	-3.8999	0.07	
18	2018-06-12 11:10:27 UTC-5	1806004-04	KF SAMPLE				0.04	0.42	-0.7999	0.09	
19	2018-06-12 11:02:50 UTC-5	B18F011-DUP2	KF SAMPLE				1.16	11.58	-0.7991	1.92	4
20	2018-06-12 10:57:30 UTC-5	1805008-02	KF SAMPLE				1.05	10.52	-0.7945	1.74	
21	2018-06-12 10:51:16 UTC-5	Hydranal 10.0 (B18F011-SRM4)	KF SAMPLE				1.04	10.44	-1.0028	2.18	
22	2018-06-12 10:47:37 UTC-5	Hydranal 10.0 (B18F011-SRM3)	KF SAMPLE				1.03	10.30	-1.0454	2.24	
23	2018-06-12 10:44:22 UTC-5	Hydranal 10.0 (B18F011-SRM2)	KF SAMPLE				1.04	10.39	-1.0347	2.23	
24	2018-06-12 10:33:39 UTC-5	Hydranal 10.0 (B18F011-SRM1)	KF SAMPLE				1.04	10.41	-1.0436	2.26	
25	2018-06-12 10:27:18 UTC-5	Water	KF TITER	4.85	4.81	0.72			-0.0521	10.74	
26	2018-06-12 10:23:37 UTC-5	Water	KF TITER	4.79	4.79	0.00			-0.0529	11.04	
27	2018-06-12 10:20:16 UTC-5	Water	KF TITER	4.79	4.79	invalid			-0.0361	7.54	






REVIEWED

By Colin Breslin at 12:19 pm, Jun 13, 2018

# Summary of Comments on Titel

---

Page: 1

	Number: 1	Author: cbreslin	Subject: Inserted Text	Date: 6/13/2018 11:34:31 AM
	repeated with larger sample size.			
	Number: 2	Author: cbreslin	Subject: Inserted Text	Date: 6/13/2018 11:33:40 AM
	over range. repeated with smaller sample size.			
	Number: 3	Author: cbreslin	Subject: Inserted Text	Date: 6/13/2018 11:32:55 AM
	repeated with larger sample size.			
	Number: 4	Author: cbreslin	Subject: Inserted Text	Date: 6/13/2018 11:32:16 AM
	repeated with larger sample size.			
	Number: 5	Author: cbreslin	Subject: Reviewed	Date: 6/13/2018 12:19:51 PM

## Flash Point Calculations

Instructions for Use:

- 1) Enter the unique sample ID under "LIMS ID".
- 2) Enter the observed flash point temp (°F) under "Flash Point Temp Recorded (°F)".
- 3) Enter the observed barometric pressure under the appropriate column to match the units recorded.
- 4) The spreadsheet automatically calculates the barometrically corrected flash point based on the units of pressure used.
- 4) The final calculated result is given under "Calculated Flash Point (°F)".
- 5) Follow CRL SOP AIG048A for reporting the final results.

Date of Analysis: 6/6/2018

Analyst: CB

Work Order(s): 1805008

LIMS Batch ID: B18F004

Count	LIMS ID	Flash Point Temp Recorded (°F)	Barometric Pressure Recorded (mm Hg)	Barometric Pressure Recorded (mb)	Barometric Pressure Recorded (in Hg)	Calculated Flash Point (°F)
1	B18F004-SRM1	79		990.0		80.05
2	B18F004-SRM2	80		990.0		81.05
	<b>B18F004-SRM3 (p-Xylene Phillips 66), average =</b>					<b>80.55</b>
3	1805008-02	145		991.0		146.00
4	1805008-02	146		990.0		147.05
	<b>1805008-02, average =</b>					<b>146.52</b>
5	B18F004-DUP1	147		988.0		148.14
6	B18F004-DUP1	147		988.0		148.14
	<b>B18F004-DUP1, average =</b>					<b>148.14</b>

## Flash Point Calculations

Instructions for Use:

- 1) Enter the unique sample ID under "LIMS ID".
- 2) Enter the observed flash point temp (°F) under "Flash Point Temp Recorded (°F)".
- 3) Enter the observed barometric pressure under the appropriate column to match the units recorded.
- 4) The spreadsheet automatically calculates the barometrically corrected flash point based on the units of pressure used.
- 4) The final calculated result is given under "Calculated Flash Point (°F)".
- 5) Follow CRL SOP AIG048A for reporting the final results.

Date of Analysis: 6/7/2018

Analyst: CB

Work Order(s): 1805008

LIMS Batch ID: B18F004

Count	LIMS ID	Flash Point Temp Recorded (°F)	Barometric Pressure Recorded (mm Hg)	Barometric Pressure Recorded (mb)	Barometric Pressure Recorded (in Hg)	Calculated Flash Point (°F)
1	B18F004-SRM4	80		990.0		81.05
2	B18F004-SRM5	80		990.0		81.05
	<b>B18F004-SRM6 (p-Xylene Phillips 66), average =</b>					<b>81.05</b>
3						45.60
4						45.60

## Flash Point Calculations

Instructions for Use:

- 1) Enter the unique sample ID under "LIMS ID".
- 2) Enter the observed flash point temp (°F) under "Flash Point Temp Recorded (°F)".
- 3) Enter the observed barometric pressure under the appropriate column to match the units recorded.
- 4) The spreadsheet automatically calculates the barometrically corrected flash point based on the units of pressure used.
- 4) The final calculated result is given under "Calculated Flash Point (°F)".
- 5) Follow CRL SOP AIG048A for reporting the final results.

Date of Analysis: 6/8/2018

Analyst: CB

Work Order(s): 1805008, 1806004

LIMS Batch ID: B18F004

Count	LIMS ID	Flash Point Temp Recorded (°F)	Barometric Pressure Recorded (mm Hg)	Barometric Pressure Recorded (mb)	Barometric Pressure Recorded (in Hg)	Calculated Flash Point (°F)
1	B18F004-SRM7	80		996.0		80.78
2	B18F004-SRM8	79		997.0		79.73
	<b>B18F004-SRM9 (p-Xylene Phillips 66), average =</b>					<b>80.25</b>
3	1806004-01	135		996.0		135.78
4	1806004-01	135		997.0		135.73
	<b>1806004-01, average =</b>					<b>135.75</b>
5	B18F004-DUP2	135		997.0		135.73
6	B18F004-DUP2	135		997.0		135.73
	<b>B18F004-DUP2, average =</b>					<b>135.73</b>
7	1806004-02	135		995.0		135.82
8	1806004-02	135		995.0		135.82
	<b>1806004-02, average =</b>					<b>135.82</b>

# Precision and Accuracy/IDOC Study

**Analyst:** Colin Breslin  
**Analysis Date:** 6/12/2018  
**SOP No. and Version:** AIG015A V4  
**SOP Name:** Water Content by Karl Fischer Titration  
**Current Limits:** 97 - 107 %

*\*To use template for calculations,  
 save the sheet as another name and  
 enter the information needed.  
 Change units and significant digits  
 based on SOP.\**

Sample LIMS ID	Spiked Conc. (%)	Measured Conc. (%)	% Recovery
B18F011-SRM1	1.00	1.04	104.0
B18F011-SRM2	1.00	1.04	104.0
B18F011-SRM3	1.00	1.03	103.0
B18F011-SRM4	1.00	1.04	104.0

CALCULATIONS								
n	$\bar{X}$	s	2s	3s	$\bar{X} - 2s$	$\bar{X} + 2s$	$\bar{X} - 3s$	$\bar{X} + 3s$
4	1.0375	0.0050	0.0100	0.0150	1.0275	1.0475	1.0225	1.0525
4	103.8	0.500	1.0	1.5	102.8	104.8	102.3	105.3

Where:

n = Number of replicate analyses  
 s = Sample (n-1) standard deviation of the replicate analyses  
 2s = Warning limits for precision control chart  
 3s = Action limits for precision control chart  
 $\bar{X}$  = Mean of the replicate analyses  
 $\bar{X} \pm 2s$  = Warning limits for accuracy control chart  
 $\bar{X} \pm 3s$  = Action limits for accuracy control chart