### **ARTICLE**

# **Evaluation of Dietary Supplement Contamination** by Xenobiotic and Essential Elements Using Microwave-Enhanced Sample Digestion and Inductively Coupled Plasma-Mass Spectrometry

Gregory M. Zinn<sup>1</sup>, G. M. Mizanur Rahman<sup>2</sup>, Scott Faber<sup>3</sup>, Mesay Mulugeta Wolle<sup>1</sup>, Matt Pamuku<sup>2</sup>, & H. M. Skip Kingston<sup>1</sup>

<sup>1</sup>Department of Chemistry and Biochemistry, Duquesne University, Pittsburgh, Pennsylvania, USA, <sup>2</sup>Applied Isotope Technologies, Pittsburgh, Pennsylvania, USA, <sup>3</sup>Department of Medicine, The Children's Institute, Pittsburgh, Pennsylvania, USA

**ABSTRACT.** Dietary supplements were analyzed by evaluating the elemental content in six widely consumed products manufactured by four well-known companies. The elements included the neurotoxic and carcinogenic elements cadmium, mercury, aluminum, lead, arsenic, and antimony, as well as the essential elements zinc, selenium, chromium, iron, and copper, which were often not listed as ingredients on the product labels. Contamination from either xenobiotic or essential elements was found in all samples analyzed. The samples were prepared using US Environmental Protection Agency (EPA) Method 3052, microwave-enhanced digestion. The resulting digests were analyzed by Inductively Coupled Plasma-Mass Spectrometry based on EPA Method 6020B. The analytical protocols were validated by analyzing a multivitamin standard reference material, the National Institute of Standards and Technology Standard Reference Material 3280. The application of EPA standard methods demonstrated their utility in making accurate and precise measurements in complex matrices with multiple ingredients and excipients. In the future, the use of these methods could provide a uniform quality assurance protocol that can be implemented along with other industry guidelines to improve the production of dietary supplements.

dietary supplements, elemental contamination, EPA Method 3052, EPA Method 6020B, inductively coupled plasma-mass spectrometry, microwaveenhanced digestion

#### INTRODUCTION

Dietary supplements are the products containing vitamins, minerals, amino acids, plant materials, and/or other nutritional additives consumed to complement the human diet (Garcia-Rico et al., 2007). They are purchased by the public to resolve

Address correspondence to: Dr. Gregory M. Zinn, Department of Chemistry and Biochemistry, Duquesne University, 600 Forbes Avenue, 308 Mellon Hall, Pittsburgh, PA 15282, USA. (Email: zinng@duq.edu).





nutrient deficiencies, prevent birth defects, treat diseases and symptoms, and im-

athletic performance (Driscoll et al., 2010; Maughan, 2005; Raman et al., 2004). Supplements are widely available in retail stores, shopping malls, gyms, gas stations, convenience stores, and on the internet (Taylor, 2004). They continue to garner attention from media exposure from studies, articles, and advertisements (Radimer et al., 2004). The reasons for the increased popularity of dietary supplements include hope for the conservation of intellectual capacity, increased energy, stress relief, mitigation of high healthcare costs, and a perception that they contribute to a healthier lifestyle (Raman et al., 2004). Physicians routinely incorporate dietary supplements as part of their patient advice and clinical treatment. It was estimated that globally consumers had spent at least US\$ 60 billion on dietary supplements in 2006 (Crowley & FitzGerald, 2006; Geyer et al., 2008), with a more recent estimate suggesting that between US\$ 21 and 25 billion were spent annually in the United States alone (McCormick, 2010). Based on a National Health Interview Survey, over half the US adult population includes supplements in their daily diets (Cohen, 2009). The wide and growing consumption of dietary supplements makes it imperative for public health that existing quality control and assurance methods, some of which are decades old, are updated and streamlined utilizing advanced measurement tools that produce consistent, reliable, and quality data. Minimizing and avoiding unsafe levels of toxins and labeling errors require the quality control process to include the entire supply chain.

The expectations of the general public for dietary supplements are the same as they are for food: contribute nutrients and avoid any constituent that will have a negative health impact. Many reported instances of tainted products suggest that the public's expectations are not being met. It has been estimated that 50,000 unreported contamination incidences occur annually in dietary supplements in the United States (Cohen, 2009). These contaminants include heavy metals, bacteria, toxic plant matter, prescription medications, controlled substances, rejected drugs, and other compounds (Cohen, 2009).

The toxic effects of xenobiotic metals on human health have been well documented in the literature. These elements, including lead (Pb), mercury (Hg), cadmium (Cd), and arsenic (As), cause a wide variety of toxicological effects and symptoms affecting all organ systems, including the central nervous system (Neustadt & Pieczenik, 2007). Lead has irreversible, detrimental effects that interfere with neuronal signaling in the central nervous system, heme synthesis, and calcium (Ca) metabolism and function (Bourgoin et al., 1993; Sargent, 1994). Mercury exposure can cause kidney damage and neurological disorders (Chen et al., 2002). Cadmium results in reproductive disruption, kidney failure, and weakening of the bones (Garcia-Rico et al., 2007). Arsenic has been connected to hypertension, hypotension, cardiomyopathy, cardiac disease, and peripheral neuropathy (Jomova et al., 2011). Aluminum (Al) has been associated with Alzheimer's disease (Kawahara & Kato-Negishi, 2011).

In parallel with the negative effects of toxic heavy metals, there is concern about improper labeling of dietary supplements that omits the presence of essential elements, such as zinc (Zn), selenium (Se), chromium (Cr), iron (Fe), and copper (Cu). Although they may not be listed as supplement ingredients, their presence



can cause adverse health consequences, including the disruption of nutritional regimens. Adulteration related to these essential elements must be considered when evaluating dietary supplements for elemental contamination in addition to xenobiotic elements. Essential elements are common ingredients of dietary supplements because they are needed to maintain crucial metabolism in the body (Raghunath et al., 2006). Human metabolism requires many different metal coenzymes for the optimal functionality of numerous proteins. Although they are necessary, these elements can have negative consequences if present in excess. Zinc is required for the normal maturation and optimal performance of the immune system (Shankar & Prasad, 1998); however, patients that exceed the daily recommended allowance for Zn could experience leukopenia, neutropenia, anemia, sideroblastic anemia, and hypocupremia (Driscoll et al., 2010; Igic et al., 2002; Porea et al. 2000; Salzman et al. 2002). While Se exists naturally and is needed for proper cellular functioning (MacFarquhar et al., 2010), it is toxic at high levels, resulting in vomiting, nausea, pulmonary edema, renal and liver damage, and death (Fan & Kizer, 1990; Nuttall, 2006). Chromium performs a key function in the metabolism of carbohydrates and fats, and increases cellular response to insulin (Marrero et al., 2013); however, an overabundance can cause kidney failure, improper liver functioning, and gastrointestinal bleeding (MacFarquhar et al., 2010). Iron is well known for promoting oxygen transport and reducing oxygen shortages in cells, thereby preventing exhaustion and a compromised immune system (Marrero et al., 2013). Excess Fe can negatively affect the metabolism of other closely related elements such as Zn, Cu, and Ca, and may result in oxidative stress, digestive disturbances, and vomiting (Fairweather-Tait et al., 2011). Copper plays a role in many essential reactions needed for the survival of mammals (Linder & Hazegh-Azam, 1996), but Cu surplus can cause liver damage (Lopez de Romana et al., 2011). It has been hypothesized that the consumption of inorganic Cu in dietary supplements might contribute to the development of Alzheimer's disease (Brewer, 2009, 2011).

In spite of the known effects of xenobiotics, the regulation of dietary supplements has become less stringent in the United States after the passing of the Dietary Supplement Health and Education Act (DSHEA) in 1994 (Cohen, 2009) as a revision to the Federal Food, Drug, and Cosmetic Act (Saldanha et al., 2004). The testing requirements for the production and sale of dietary supplements were significantly relaxed in the DSHEA (Maughan, 2005). Dietary supplements were previously classified as being "food additives," which had mandated evidence of their safety before being released to consumers (Cohen, 2009). The DSHEA created a special category of food called "dietary supplements" (Saldanha et al., 2004), with provisions that were not applicable to conventional and functional foods (Taylor, 2004). Since dietary supplements do not diagnose, prevent, or cure disease according to the US Congress, they are not required to have the same regulatory standards as pharmaceutical drugs (Maughan et al., 2011). Subsequently, there is no US Food and Drug Administration (FDA)-required testing of dietary supplements for safety and efficacy prior to market release, and approval is not needed before the manufacture and distribution of these products (Ashar et al., 2007). Under the current guidelines, the FDA will only recommend testing after a dietary supplement has been determined to be unsafe (Woo, 2007). This relaxed regulatory requirement is



very different from the treatment of prescription drugs, which require FDA review of safety and effectiveness before reaching the market (Ashar et al., 2007). Without providing specific guidelines of dietary supplement production, the DSHEA established Good Manufacturing Practices (GMPs) to maintain the quality of dietary supplements available to the general public (Ashar et al., 2007). Extensive control over the entire manufacturing process is required for GMP compliance, including starting ingredients, finished goods, and quality control (Tumir et al., 2010). The DSHEA does not specify standard analytical methods, such as those from the US Environmental Protection Agency (EPA), which can be adopted and implemented for quality assurance as a way to improve public safety and attain GMP compliance. The US EPA is a globally recognized and respected repository of workgroup-tested and optimized methods, which may be adopted as templates and optimized as standard quality assurance and testing protocols.

Different government agencies have limits for exposure to certain elements referred to as maximum contaminant levels (MCLs), which are shown for each agency in Table 1. The United States Pharmacopeia (USP) has set elemental maximum values in dietary supplements for only inorganic As, Cd, Pb, total Hg, and methylmercury (United States Pharmacopeia, 2012b). This amount of regulation is much less stringent than what the USP requires for prescription drugs, where 15 elements are listed with concentration thresholds (United States Pharmacopeia, 2012a). The FDA does not have specific elemental limits for dietary supplements but rather focuses on GMPs to prevent the adulteration of dietary supplements from heavy metals, pesticides, and other harmful contaminants, and to ensure the presence of only the listed ingredients (Mindak et al., 2008). California Proposition 65 regulations provide a list of chemicals and compounds which are known to cause cancer or reproductive toxicity. Daily limits (Office of Environmental Health Hazard Assessment, 2012) and MCLs (Office of Environmental Health Hazard Assessment, 2012) are provided for some elements, but there is no specific information for dietary supplements. The EPA provides MCLs in drinking water, but the amount of drinking water consumed per day is typically much greater than the amount of dietary supplements ingested (US Environmental Protection Agency (EPA)). Recent papers investigating elemental contamination in dietary supplements have compiled tables with the advocated amounts of selected elements to be consumed each day (Avula et al., 2010, 2011). Commonly, these limits are calculated based on the average adult male; often they do not take into account women, children, and developing fetuses that can be more susceptible to the effects of these potentially harmful elements and chemicals. It is hard to assess the quality of dietary supplements based on the current guidelines provided by these regulatory agencies. However, the USP-published methods do not match the level of completeness and utility of the EPA methods as many USP methods are either too brief, inconsistent, or too old to be practical or useful today to achieve GMP compliance.

Numerous studies have demonstrated the problem of harmful elemental contamination in dietary supplements. Nearly 20 years ago, Pb was discovered in calcium supplements (Bourgoin et al., 1993), while a more recent study determined that two-thirds of calcium supplements contained Pb concentrations that exceeded the 1999 California limits (Scelfo & Flegal, 2000). Another study found wide concentration ranges for As, Cd, Hg, and Pb in dietary supplements (Dolan et al.,



 TABLE 1. List of Inorganic Contaminants and Their MCLs by Regulatory Agency

Element	USP Dietary Supplement Limit ( $\mu$ g g <sup>-1</sup> ) (United States Pharmacopeia 2012b)	USP Drug Product Limit ( $\mu$ g g <sup>-1</sup> ) (United States Pharmacopeia 2012a)	Prop 65 $NSRL^a$ or $MADL^b$ ( $\mu g day^{-1}$ ) (Office of Environmental Health Hazard Assessment)	Prop 65 MCL (mg L <sup>-1</sup> ) (Office of Environmental Health Hazard Assessment)	EPA Drinking Water MCL (mg L <sup>-1</sup> ) (United States Environmental Protection Agency (EPA))
Aluminum	No limit	No limit	No limit	-	0.05 to 0.2
Antimony	No limit	No limit	No limit	0.006	0.006
Arsenic	1.5	0.15	0.06 (inhalation) 10 (except inhalation)	0.01	0.010
Barinm	No limit	No limit	No limit	-	2
Beryllium	No limit	No limit	0.1	0.004	0.004
Cadmium	0.5	2.5	0.05 (inhalation) 4.1 (oral)	0.005	0.005
Chromium (total)	No limit	No limit	No limit	No limit	0.100
Copper	No limit	100	No limit	1.3	1.3
Iridium	No limit	10	No limit	No limit	No limit
Iron	No limit	No limit	No limit	No limit	0.3
Lead	1.0	0.5	15 (oral) 0.5	0.015	0.015
Manganese	No limit	No limit	No limit	No limit	0.05
Mercury	1.5	1.5	No limit	0.002	0.002
(Inorganic)					
Molybdenum	No limit	10	No limit	No limit	No limit
Nickel	No limit	50	No limit	0.1	No limit
Osmium	No limit	10	No limit	No limit	No limit
Palladium	No limit	10	No limit	No limit	No limit
Platinum	No limit	10	No limit	No limit	No limit
Rhodium	No limit	10	No limit	No limit	No limit
Ruthenium	No limit	10	No limit	No limit	No limit
Selenium	No limit	No limit	No limit	No limit	0.05
Silver	No limit	No limit	No limit	No limit	0.10
Thallium	No limit	No limit	No limit	0.002	0.002
Uranium	No limit	No limit	No limit	25	0:030
Vanadium	No limit	10	No limit	No limit	No limit
Zinc	No limit	No limit	No limit	No limit	5

aNSRL: No significant risk level; b MADL: maximum allowable dose level.



2003). Studies from multiple countries noted xenobiotic and essential elemental contamination in natural drugs in China (Chuang et al., 2000), in herbal remedies in Nigeria (Obi et al., 2006), and in dietary supplements in Mexico (Garcia-Rico et al., 2007) and the United States (Avula et al., 2010, 2011).

Other researchers have investigated contamination from essential elements in dietary supplements. A recent study discussed widespread Se toxicity resulting from elevated Se quantities in a liquid dietary supplement that also contained a high concentration of Cr (MacFarguhar et al., 2010). Additional case studies have documented the toxicity caused by Se in dietary supplements (Sutter et al., 2008; Aldosary et al., 2012).

The project described in this paper started with an announcement from a US dietary supplement supplier that one of their supplements had been contaminated. This study expanded on the original discovery by evaluating elemental content in six common types of dietary supplements (zinc, essential fatty acids, calcium, magnesium, multivitamins, and probiotics) from four different manufacturers. The present study applied additional quality control parameters to the elemental analysis of dietary supplements through the use of the US Environmental Protection Agency (EPA) (1996) Method 3052 for sample preparation and the US Environmental Protection Agency (EPA) (2013) Method 6020B for analysis, with National Institute of Standards and Technology (NIST) Multivitamin/Multielement Tablets, Standard Reference Material (SRM) 3280, used for validation.

### MATERIALS AND METHODS

# Materials

### Chemicals and Standards

All reagents were of analytical grade. ARISTAR ULTRA nitric acid (HNO<sub>3</sub>), ARISTAR ULTRA hydrochloric acid (HCl), ARISTAR ULTRA water, borosilicate glass vessels, and 50-ml polypropylene centrifuge tubes were procured from VWR (Radnor, PA, USA). Multi-element and single element standards of Hg, iodine (I), lithium (Li), and yttrium (Y) were bought from Inorganic Ventures (Christiansburg, VA, USA). The tuning mix for the Inductively Coupled Plasma-Mass Spectrometer (ICP-MS) comprising 0.001- $\mu g$  g<sup>-1</sup> lithium, cobalt, yttrium, cerium, and tellurium in 2% (v:v) HNO<sub>3</sub> was purchased from Agilent Technologies (Santa Clara, CA, USA). The internal standard solution containing 100 µg g<sup>-1</sup> of <sup>6</sup>lithium, scandium, germanium, rhodium, indium, terbium, lutetium, and bismuth was also obtained from Agilent Technologies. Ultra high purity grade (99.999%) liquid argon along with hydrogen and helium gases was procured from Airgas (Radnor, PA, USA). SRM 3280 was bought from NIST (Gaithersburg, MD, USA).

## Dietary Supplement Samples

Six types of dietary supplements were purchased from four manufacturers in commercially available form (Table 2). From Company 1, zinc, calcium, and magnesium supplements were obtained. Zinc, essential fatty acids, calcium, magnesium,



and multivitamin supplements were tested from Company 2. All six supplement types were analyzed from Company 3 and 4.

### Instrumentation

Sample preparation was accomplished by applying microwave-enhanced chemistry using a Milestone (Shelton, CT, USA) UltraWAVE laboratory microwave system equipped with temperature and pressure feedback control. This device accurately monitored and controlled temperature to within  $\pm 2^{\circ}$ C of the specified values, automatically adjusting the microwave field output power to achieve the preset temperatures. The UltraWAVE allowed for extremely high temperature and pressure conditions of up to 300°C and 199 bar, respectively, to be maintained for extended periods during the digestion of the samples. This capability allowed for different types of samples to be digested simultaneously since all samples were held under the same conditions. In this study, borosilicate glass vessels were used for the digestion of 15 samples per batch.

Mass spectrometric analysis was performed using an Agilent 7700 ICP-MS with a micromist nebulizer, quartz spray chamber, an octopole reaction system, and a quadrupole mass analyzer. This instrument has a collision/reaction cell that uses hydrogen or helium gas for the elimination of polyatomic interferences. All samples and the internal standard were introduced into the ICP-MS using the peristaltic pump at a speed of 0.1 revolutions per second (rps) via the micromist nebulizer of the spray chamber. Before analysis, the samples were placed in a CETAC (Omaha, NE, USA) ASX-520 autosampler housed inside an ENC-500 anti-contamination enclosure. The ICP-MS and autosampler were located in an International Organization for Standardization Class 6 cleanroom laboratory. The optimized parameters of the ICP-MS and autosampler are shown in Table 3.

### **METHODS**

# Homogenization of the Dietary Supplements

Four different forms of dietary supplements were analyzed during this study: tablets, capsules, powders, and liquids. The procedure given in the Certificate of Analysis for NIST SRM 3280 was followed for sample homogenization. Fifteen SRM 3280 tablets were homogenized in a clean mortar and pestle and stored in a 50-ml polypropylene centrifuge tube. All of the tablet dietary supplement samples were prepared and stored in the same manner as SRM 3280. For each capsule sample, the outer casing was removed for 15 capsules and the inside contents were emptied into a 50-ml polypropylene centrifuge tube. For the powder samples, a quantity similar to 15 SRM 3280 tablets was poured into a 50-ml polypropylene centrifuge tube. The liquid samples were shaken immediately before being added to the microwave digestion vessels.

# Microwave-Enhanced Digestion of the Dietary Supplements by EPA Method 3052

The homogenized samples were digested according to the EPA Method 3052 (Kingston et al., 1997; US Environmental Protection Agency (EPA), 1996) to determine their total elemental content. In each glass vessel, 0.25 g of sample was



The Dietary Supplements Analyzed in the Study TABLE 2.

	odppicing	incording approximation of the control of the contr	
Manufacturer	Supplement	Nutritional Ingredient Type	Other Ingredients
Company 1	Zinc Calcium Magnesium	Mineral Mineral Mineral	Cellulose, water Ascorbyl palmitate, cellulose, water Ascorbyl palmitate cellulose water
Company 2	Zinc Zinc Essential fatty acid	Mineral  Mineral  Omega-3 fatty acids (docosahexaenoic acid, eicosapentaenoic acid)	Plant cellulose, stearic acid, rice flour, carob Cholesterol, vitamins A, D <sub>3</sub> , E, natural lemon flavor
	Calcium	Mineral	Fructose, sorbitol, stearic acid, magnesium stearate, vanilla, silicon dioxide
	Magnesium Multivitamin	Vitamins and minerals (vitamins A, C, D <sub>3</sub> , E, thiamin, riboflavin, niacin, vitamin B6, folate, vitamin B12, biotin, pantothenic acid, calcium, iodine, magnesium, zinc, selenium, copper,	Gelatin capsule Microcrystalline cellulose, dicalcium phosphate, croscarmellose sodium, stearic acid, magnesium stearate, ethylcellulose, hydroxypropyl cellulose
Company 3	Zinc Essential fatty	Mineral Fatty acid (marine oil concentrate,	Vitamins C and B6, cellulose, vegetarian capsule, L-leucine Natural mixed tocopherol vitamin E from soy, natural lemon
	acid Calcium	docosahexaenoic acid, eicosapentaenoic acid) Mineral	oil, rosemary extract Phosphorus, hydroxypropyl methylcellulose, water, L-leucine, cellulose
	Magnesium Multivitamin	Mineral Vitamins and minerals (vitamins A, C, D3, E, K1, thiamine, riboflavin, niacin, vitamin B6, folate, vitamin B12, biotin, pantothenic acid, calcium, iodine, magnesium, zinc, selenium, manganese,	Hydroxypropyl methylcellulose, water, L-leucine Natural citrus flavor
	Probiotic	vanadium, molyboerdm, potassium, boron, vanadium, choline, inositol, citrus bioflavonoids) Probiotic (probiotic blend, <i>lactobacillus</i> acidophilus, bifidobacterium bifidum)	Hydroxypropyl methylcellulose, water, L-leucine



Purified water, sorbitol, glycerine, citric acid, natural raspberry flavor, potassium sorbate	Cod liver oil, purified water, gelatin, glycerine	Calcium chelate, sodium starch, glycolate, magnesium stearate, vitamin D3, plant cellulose capsule	Magnesium glycinate, L-leucine, plant cellulose capsule	Dextrose, calcium carbonate, ascorbic acid, vitamin E acetate, magnesium oxide, zinc citrate, magnesium	glycinate, calcium citrate, niacinamide, d-calcium	pantothenate, manganese citrate, CoQ10, citric acid,	natural peach mango flavor, starch, stearic acid, silicon	dioxide, magnesium stearate, biotin, vitamin A acetate,	pyridoxine, hydrochloride, riboflavin, thiamine mononitrate,	glycerides of stearic and palmitic acids, mannitol, vitamin	D-3, selenium chelate, di-calcium phosphate, beta	carotene, cyanocobalamin, chromium chelate,	molybdenum chelate, sucralose, folic acid, potassium	iodide	Stabilized probiotic strains, microcrystalline cellulose, silicon	dioxide, magnesium stearate, plant cellulose capsule		
Mineral	Omega-3 fatty acid and vitamins (docosahexaenoic acid, eicosapentaenoic acid vitamin A, vitamin D)	Mineral	Mineral	Vitamins and minerals (vitamins A, C, D, E, B1, B2, niacinamide, B6, folic acid, B12, biotin,	pantothenic acid, calcium, iodine, magnesium,	zinc, selenium, manganese, chromium,	molybdenum, coenzyme Q10)								Probiotic (lactobacillus rhamnosus, lactobacillus	acidophilus, bifidobacterium bifidum/lactis,	lactobacillus casei, lactobacillus plantarum,	streptococcus therophilus)
Zinc	Essential fatty acid	Calcium	Magnesium	Multivitamin											Probiotic			
Company 4																		



TABLE 3. Optimized Operating Conditions for the Agilent 7700 ICP-MS for the Analysis of Dietary Supplement Samples According to EPA Method 6020B

Agilent 7700 ICP-MS Parameter	EPA Method 6020B Analysis Setting
Rf power	1,550 W
Spray chamber temperature	2°C
Acquisition mode	Spectrum
Peak pattern	. 3 points mass⁻¹
Sample depth	8 mm
Plasma gas	15 L min <sup>−1</sup>
Carrier gas	0.9 L min <sup>-1</sup>
Makeup gas	0.15 L min <sup>-1</sup>
H <sub>2</sub> cell gas flow	$4.5 \text{ ml min}^{-1}$
He cell gas flow	5.5 ml min <sup>-1</sup>
Gas stabilization time	30 s
Integration time	$0.1-1.0 \text{ s mass}^{-1}$
Replicates	4
Nebulizer pump speed	0.1 rps
Uptake time	30 s
Sabilization time	30 s

combined with 8-ml HNO<sub>3</sub> and 2-ml HCl. Three subsamples were prepared for each dietary supplement. The samples were digested in the microwave following a program that implemented a 10-min ramp to 180°C with a 10-min hold at 180°C. Complete digestion for all the types of dietary supplements (tablets, capsules, powders, and liquids) was achieved resulting in the production of homogeneous, transparent, light yellow solutions that did not contain any solids. After digestion, each sample was poured into a 50-ml polypropylene centrifuge tube and diluted with water to approximately 20 ml. The polypropylene centrifuge tubes were weighed before and after the addition of the sample digests. Six analytical blanks comprising HNO<sub>3</sub> and HCl were digested with the samples. The samples were stored in a cold room at 4°C and usually analyzed within three days of digestion. The microwave vessel caps were cleaned after digestion by soaking them in 1% HNO<sub>3</sub> overnight to prevent cross contamination.

# Elemental Concentration Determination by ICP-MS Using EPA Method 6020B

The total concentrations of elements in the sample digests were determined according to the EPA Method 6020B (US Environmental Protection Agency (EPA), 2013). On each day of analysis, the system was aspirated with 2% HNO<sub>3</sub> for 30 min before tuning the instrument using the Agilent tuning solution at a concentration of  $0.001 \,\mu\mathrm{g}\,\mathrm{g}^{-1}$  and determining the pulse/analog factor. Analysis on the Agilent 7700 ICP-MS was performed in three modes: hydrogen, helium, and no gas. The amount of hydrogen and helium gas was optimized to eliminate polyatomic interferences. Table 4 lists the analysis mode for each element. Between analysis modes, a 30-s delay was employed to equilibrate the presence or absence of collision gas and/or reaction gas.

Due to the wide concentration range of the analytes of interest in SRM 3280 and the dietary supplement samples, calibration standards were prepared for 60 elements at three different concentration ranges (low, medium, and high) of 0–0.025,



TABLE 4. The Analysis Information, Correlation Coefficient, Detection Limits (DLs), and Background Equivalent Concentrations (BECs) for the Xenobiotic and Essential Elements of Interest Determined by the Agilent 7700 ICP-MS

Element	M/z	Mode	r <sup>2</sup>	DL (ng $g^{-1}$ )	BEC (ng g <sup>-1</sup> )
Mg	24	He	0.9976	25.7	13.2
Al	27	$H_2$	0.9995	3.44	0.623
Ca	40	$H_2^-$	0.9983	110	93.9
V	51	He	0.9999	0.0537	0.0387
Cr	52	He	0.9999	0.961	3.76
Mn	55	He	1.0000	0.346	0.284
Fe	56	He	0.9985	18.8	15.4
Co	59	He	1.0000	0.0197	0.0438
Ni	60	He	0.9994	0.411	1.78
Cu	63	He	0.9999	0.309	0.684
Zn	66	He	0.9979	9.35	6.39
As	75	He	0.9999	0.0187	0.334
Se	78	$H_2$	1.0000	0.0178	0.0232
Мо	95	He	1.0000	0.102	0.334
Cd	111	He	1.0000	0.0477	0.0599
Sn	118	He	0.9999	0.105	0.295
Sb	121	No gas	0.9998	0.0866	0.130
Hg	201	No gas	1.0000	0.0658	0.0246
Pb	208	No gas	0.9999	0.00721	0.0752
U	238	No gas	1.0000	0.000834	0.00175

0-0.050, and 0-0.150  $\mu g$  g<sup>-1</sup>, respectively. The calibration standards were prepared by acid matching, which consisted of adding the appropriate percentage of microwave-digested reagent blank to each calibration standard.

The sample digests were diluted with water to total factors of 320, 8,000, and 800,000 for the low, medium, and high calibration ranges, respectively. These were analyzed in three separate sequences with four replicates taken from each subsample. Analyses were performed in spectrum mode using three points per mass. The integration time for each element ranged from 0.1–1.0 s per point depending on the element's ionization energy. The measured concentration value for each element in SRM 3280 and each dietary supplement sample was determined based on the value that fell within the range of one of the three calibration curves. After analysis was completed, the data were exported from the Agilent MassHunter software (version G7201A A.01.01) into Microsoft Excel for further data analysis.

# Quality Control

The current study implemented numerous quality control measures mentioned previously that included the preparation of multiple subsamples, replicate analysis, preparation of analytical blanks, acid-matched calibration standards, and SRM 3280. It also utilized many of the measures described in the EPA Method 6020B (US Environmental Protection Agency (EPA), 2013.) An Agilent internal standard solution was used at a concentration of 1  $\mu$ g g<sup>-1</sup> in 1% HNO<sub>3</sub> and 0.5% HCl, and its acceptable recovery was maintained between 80 and 120% of the calibration blank. Between each sample, the autosampler probe was washed in three solutions of 1% HNO<sub>3</sub> and 0.5% HCl for 30 s each with a nebulizer pump speed of 0.5 rps.



After every nine samples, low- and mid-point calibration checks along with an SRM 3280 assessment were performed for additional calibration verification to ensure the consistency of the results during the complete duration of the sample run.

### RESULTS

# Analytical Figures of Merit and Validation of EPA Methods

All elements displayed linearity on the six- or seven-point external calibration curves for each of the three calibration ranges (Table 4). They had correlation coefficient (r<sup>2</sup>) values greater than the 0.998 specified in the EPA Method 6020B except for magnesium (Mg, 0.9976) and Zn (0.9979) being just below the cutoff value only for the low calibration range. When rounded to three significant figures, they meet the specified requirements of the EPA Method 6020B. The four replicates for each calibration point displayed a variation of less than 5%, indicating that the obtained measurements were highly reproducible.

The elemental detection limits (DL) were defined as the average of the blank concentration for each element plus three times the standard deviation of the blank measurements for that element (Skoog et al., 2007). The DLs for all the elements in Table 4 were below  $0.001 \,\mu g \, g^{-1}$  except for Al, Ca, Fe, Mg, and Zn. The background equivalent concentrations (BECs) were calculated by multiplying the background intensity signal and the concentration of the analyte and dividing by the difference of the analyte intensity minus the intensity of the background (Thomas, 2008). The BECs were under 0.001  $\mu g g^{-1}$  for all elements other than Ca, Cr, Fe, Mg, Ni, and Zn (Table 4).

Additional confidence in the accuracy and precision of the results of the present study was achieved by using SRM 3280 to validate the sample preparation and analysis protocols. This SRM is certified for both xenobiotic and essential elements with their concentrations representing the wide range of values seen in commercially available products. Results obtained for SRM 3280 during this study are shown in Table 5 with the mean (n = 12) and 95% confidence intervals (CI). The measured values are in strong agreement with the certified and reference values, within 10% for nearly every element. The measured values overlapped the certified or reference CIs for 14 elements with the intervals for two other elements nearly intersecting, which is the same number of elements provided in two other studies (Avula et al., 2010, 2011). The recovery efficiencies ranged from 68.2 to 124% based on a comparison of the mean-listed values to the average measured quantities. Only three of the 23 certified or reference elements had either recoveries that differed by more than 20% from the mean value or that did not overlap CI values. This study provided measured values for the certified elements As, Cd, and Pb in SRM 3280 that were not included in two other studies (Avula et al., 2010, 2011). The value for As fell within the provided CI value, while Cd had a 118% recovery and Pb had a 70.7% recovery. The accuracy of all measurements in this study improved with the inclusion of the provided confidence intervals.

The precision determined by the 95% CIs was within 10% of the measured mean values for all of the certified or reference elements in SRM 3280, with over half inside 5% (Table 5). This level of precision for SRM 3280 was maintained for the



TABLE 5. Concentrations of All the Certified and Reference Elements Determined in SRM 3280 (N = 12,95% CIs)

Element	Certified Value ( $\mu$ g g <sup>-1</sup> )	Measured Value ( $\mu$ g g <sup>-1</sup> )
В	141 ± 7	121.9 ± 4.7
Mg	$67,800 \pm 4,000$	$72,600\pm2,320$
P	$75,700\pm3,200$	$80,390\pm 1,640$
CI	$53,000\pm2,300$	Not Determined
K	$53,100\pm7,000$	$65,670\pm4,370$
Ca	$110,700\pm 5,300$	$108,400\pm 8,300$
Cr	$93.7 \pm 2.7$	$96.09 \pm 2.03$
Mn	$1,440 \pm 110$	1, 627 $\pm$ 161
Fe	$12,350 \pm 910$	11, 950 $\pm$ 170
Ni	$8.43 \pm 0.30$	$5.861 \pm 0.299$
Cu	$1,400 \pm 170$	1, 591 $\pm$ 117
Zn	$10,150\pm 810$	$11,540 \pm 430$
As	$0.132 \pm 0.044$	$0.1625 \pm 0.0138$
Se	$17.42 \pm 0.45$	$17.98 \pm 0.37$
Mo	$70.7 \pm 4.5$	$84.12 \pm 1.16$
Cd	$0.08015 \pm 0.00086$	$0.09579 \pm 0.00828$
	$132.7 \pm 6.6$	Not Determined
⊃b	$0.2727 \pm 0.0024$	$0.1973 \pm 0.0077$
	Reference Value ( $\mu$ g g <sup>-1</sup> )	Measured Value ( $\mu$ g g <sup>-1</sup> )
Na	$330\pm20$	$301.7 \pm 16.7$
Si	$2,010 \pm 10$	Not Determined
Ti	$5,400 \pm 300$	Not Determined
V	8 ± 2	$9.691 \pm 0.635$
Co	$0.81\pm0.01$	$0.8572 \pm 0.0167$
Sr	$29.8 \pm 0.2$	$29.96 \pm 1.24$
Sn	$11.1 \pm 0.9$	$11.74 \pm 0.62$
Sb	$0.159 \pm 0.008$	$0.1085 \pm 0.0057$
La	$0.70\pm0.01$	$0.6151 \pm 0.0118$

remainder of the study as most results had 95% CIs that were less than 10% of their mean value. A large number of these were lower than 5% and some approached 1%. The other two studies (Avula et al., 2010, 2011), which have evaluated elemental content in botanicals and dietary supplements, supplied no standard deviation or confidence interval information for their sample measurements, listing only standard deviations for their SRM 3280 results. The exclusion of standard deviation and confidence intervals provided no information on their ability to sustain precise measurements for their samples during the remainder of their studies. This study provided 95% CIs for all detectable elemental quantities in the dietary supplement samples.

# Determination of Elemental Concentrations in Dietary Supplements Using EPA Method 6020B

The elemental analysis results from the four dietary supplement companies are shown in Tables 6–9, respectively. The results are presented by company and separated into xenobiotic and essential elements. Only the most important results will be described in the text form. Contamination from essential elements was considered to exist if that element was not listed as an ingredient on the label for that particular dietary supplement.



**TABLE 6.** Concentrations ( $\mu q q^{-1}$ ) of Xenobiotic and Essential Elements Determined in Dietary Supplements from Company 1 (N = 12,95% CIs)

	Zinc Supplement	Calcium Supplement	Magnesium Supplement
Xenobiotic Elements			
Al	<dl< td=""><td><dl< td=""><td><dl< td=""></dl<></td></dl<></td></dl<>	<dl< td=""><td><dl< td=""></dl<></td></dl<>	<dl< td=""></dl<>
As	<dl< td=""><td><dl< td=""><td><dl< td=""></dl<></td></dl<></td></dl<>	<dl< td=""><td><dl< td=""></dl<></td></dl<>	<dl< td=""></dl<>
Cd	$0.1168 \pm 0.0089$	<dl< td=""><td><dl< td=""></dl<></td></dl<>	<dl< td=""></dl<>
Hg	<dl< td=""><td><dl< td=""><td><dl< td=""></dl<></td></dl<></td></dl<>	<dl< td=""><td><dl< td=""></dl<></td></dl<>	<dl< td=""></dl<>
Pb	<dl< td=""><td><math>0.1718 \pm 0.0032</math></td><td><math>0.009087 \pm 0.001802</math></td></dl<>	$0.1718 \pm 0.0032$	$0.009087 \pm 0.001802$
Sb	<dl< td=""><td><dl< td=""><td><dl< td=""></dl<></td></dl<></td></dl<>	<dl< td=""><td><dl< td=""></dl<></td></dl<>	<dl< td=""></dl<>
Sn	$0.02260 \pm 0.00895$	$0.7526 \pm 0.0355$	$0.02330 \pm 0.00553$
U	<dl< td=""><td><math>0.0009512 \pm 0.0001482</math></td><td><math>0.004564 \pm 0.000150</math></td></dl<>	$0.0009512 \pm 0.0001482$	$0.004564 \pm 0.000150$
V	<dl< td=""><td><math display="block">0.02266 \pm 0.00509</math></td><td><math display="block">2.129 \pm 0.050</math></td></dl<>	$0.02266 \pm 0.00509$	$2.129 \pm 0.050$
Essential Elements			
Ca	<dl< td=""><td>257,500 <math>\pm</math> 8,600</td><td><math>208.3 \pm 14.3</math></td></dl<>	257,500 $\pm$ 8,600	$208.3 \pm 14.3$
Co	<dl< td=""><td><dl< td=""><td><math>0.05117 \pm 0.00401</math></td></dl<></td></dl<>	<dl< td=""><td><math>0.05117 \pm 0.00401</math></td></dl<>	$0.05117 \pm 0.00401$
Cr	<dl< td=""><td><dl< td=""><td><math>2.430 \pm 0.119</math></td></dl<></td></dl<>	<dl< td=""><td><math>2.430 \pm 0.119</math></td></dl<>	$2.430 \pm 0.119$
Cu	<dl< td=""><td><math>0.3349 \pm 0.0236</math></td><td><math>0.2421 \pm 0.0130</math></td></dl<>	$0.3349 \pm 0.0236$	$0.2421 \pm 0.0130$
Fe	<dl< td=""><td><math>7.974 \pm 1.439</math></td><td><math>15.81 \pm 1.63</math></td></dl<>	$7.974 \pm 1.439$	$15.81 \pm 1.63$
Mg	<dl< td=""><td><math>5,\!578 \pm 807</math></td><td>138,900 <math>\pm</math> 4,800</td></dl<>	$5,\!578 \pm 807$	138,900 $\pm$ 4,800
Mn	<dl< td=""><td><math>0.3570 \pm 0.0306</math></td><td><math>3.192 \pm 0.133</math></td></dl<>	$0.3570 \pm 0.0306$	$3.192 \pm 0.133$
Мо	<dl< td=""><td><math>0.05454 \pm 0.01049</math></td><td><math>0.5963 \pm 0.0210</math></td></dl<>	$0.05454 \pm 0.01049$	$0.5963 \pm 0.0210$
Ni	<dl< td=""><td><dl< td=""><td><math>1.774 \pm 0.073</math></td></dl<></td></dl<>	<dl< td=""><td><math>1.774 \pm 0.073</math></td></dl<>	$1.774 \pm 0.073$
Se	<dl< td=""><td><math>0.04391 \pm 0.00265</math></td><td><math>0.02158 \pm 0.00451</math></td></dl<>	$0.04391 \pm 0.00265$	$0.02158 \pm 0.00451$
Zn	$\textbf{94,200} \pm \textbf{1,100}$	$83.28 \pm 2.81$	<dl< td=""></dl<>

<sup>&</sup>lt;DL: Below detection limit.</p>

Elements in bold are not considered elemental contamination due to their inclusion in the ingredient list.

### Xenobiotic Elements

The supplements from Company 1 were found to have less contamination from xenobiotic elements compared with the other companies (Table 6). The calcium supplement contained tin (Sn) at a concentration of over 0.7  $\mu$ g g<sup>-1</sup> and Pb above  $0.17~\mu \mathrm{g}~\mathrm{g}^{-1}$ . The magnesium supplement had a similar frequency of contamination with vanadium (V) found at more than  $2 \mu g g^{-1}$  along with Sn and Pb. Most notably, the zinc supplement had 0.117  $\mu$ g g<sup>-1</sup> of Cd.

Similar results were found in the dietary supplements analyzed from Company 2 with a high frequency of contamination by xenobiotic elements (Table 7). The concentration of Al in the magnesium, zinc, and multivitamin supplements was determined to be more than 10  $\mu g$  g<sup>-1</sup>. The magnesium supplement also contained more than 7.8  $\mu$ g g<sup>-1</sup> of V. Except for the essential fatty acids supplement, Pb contamination was observed in all supplements, with the zinc supplement having a concentration of 0.738  $\mu$ g g<sup>-1</sup>, the highest of the dietary supplements tested in this study. The multivitamin supplement had the highest amount of contamination from xenobiotic elements that included As, Cd, antimony (Sb), Sn, U, and V. The zinc supplement contained Cd, Sb, Sn, and V. The calcium supplement was sullied with V and As. The essential fatty acids supplement was not contaminated with any xenobiotic elements.

Among the products tested, Company 3 had the highest frequency of xenobiotic contamination in its products (Table 8). The calcium, magnesium, zinc, and multivitamin supplements contained high amounts of Al with the mean concentra-



**TABLE 7.** Concentrations ( $\mu$ g g<sup>-1</sup>) of Xenobiotic and Essential Elements Determined in Dietary Supplements from Company 2 (N = 12,95% Cls)

	Zinc Supplement	Essential Fatty Acids Supplement	Calcium Supplement	Magnesium Supplement	Multivitamin Supplement
Xenobiotic Elements					
A	$10.75 \pm 2.92$	<dl< td=""><td><dl< td=""><td><math>61.07 \pm 15.91</math></td><td><math>12.20 \pm 3.00</math></td></dl<></td></dl<>	<dl< td=""><td><math>61.07 \pm 15.91</math></td><td><math>12.20 \pm 3.00</math></td></dl<>	$61.07 \pm 15.91$	$12.20 \pm 3.00$
As	<dl< td=""><td><dl< td=""><td><math>0.01535 \pm 0.00353</math></td><td><math>0.04374 \pm 0.00499</math></td><td><math>0.1484 \pm 0.0125</math></td></dl<></td></dl<>	<dl< td=""><td><math>0.01535 \pm 0.00353</math></td><td><math>0.04374 \pm 0.00499</math></td><td><math>0.1484 \pm 0.0125</math></td></dl<>	$0.01535 \pm 0.00353$	$0.04374 \pm 0.00499$	$0.1484 \pm 0.0125$
PO	$0.07024 \pm 0.01027$	<dl< td=""><td><dl< td=""><td><dl< td=""><td><math>0.07948 \pm 0.00890</math></td></dl<></td></dl<></td></dl<>	<dl< td=""><td><dl< td=""><td><math>0.07948 \pm 0.00890</math></td></dl<></td></dl<>	<dl< td=""><td><math>0.07948 \pm 0.00890</math></td></dl<>	$0.07948 \pm 0.00890$
Нg	<dl< td=""><td><dl< td=""><td><dl< td=""><td>&lt; DL</td><td>^DL</td></dl<></td></dl<></td></dl<>	<dl< td=""><td><dl< td=""><td>&lt; DL</td><td>^DL</td></dl<></td></dl<>	<dl< td=""><td>&lt; DL</td><td>^DL</td></dl<>	< DL	^DL
Pb	$0.7384 \pm 0.0193$	<dl< td=""><td><math display="block">0.03124 \pm 0.00598</math></td><td><math>0.05145 \pm 0.00232</math></td><td><math>0.2479 \pm 0.0097</math></td></dl<>	$0.03124 \pm 0.00598$	$0.05145 \pm 0.00232$	$0.2479 \pm 0.0097$
Sb	$0.01004 \pm 0.00166$	<dl< td=""><td><dl< td=""><td>&lt;<b>D</b>F</td><td><math>0.05058 \pm 0.00363</math></td></dl<></td></dl<>	<dl< td=""><td>&lt;<b>D</b>F</td><td><math>0.05058 \pm 0.00363</math></td></dl<>	< <b>D</b> F	$0.05058 \pm 0.00363$
Sn	$0.04709 \pm 0.00619$	<dl< td=""><td><dl< td=""><td>&lt;<u>0</u>F</td><td><math>0.09608 \pm 0.01440</math></td></dl<></td></dl<>	<dl< td=""><td>&lt;<u>0</u>F</td><td><math>0.09608 \pm 0.01440</math></td></dl<>	< <u>0</u> F	$0.09608 \pm 0.01440$
n	$0.003061 \pm 0.000240$	<dl< td=""><td><math>0.004006 \pm 0.000166</math></td><td><math>0.007497 \pm 0.000349</math></td><td><math>0.3710 \pm 0.0159</math></td></dl<>	$0.004006 \pm 0.000166$	$0.007497 \pm 0.000349$	$0.3710 \pm 0.0159$
>	$0.03149 \pm 0.00577$	<dl< td=""><td><math>0.05367 \pm 0.00458</math></td><td><math>7.849 \pm 0.114</math></td><td><math display="block">0.7924 \pm 0.0327</math></td></dl<>	$0.05367 \pm 0.00458$	$7.849 \pm 0.114$	$0.7924 \pm 0.0327$
Essential Elements					
Ca	$547.3 \pm 76.6$	<dl< td=""><td><math>85,970 \pm 1,500</math></td><td><math>815.1 \pm 56.7</math></td><td><math display="block">{\bf 43,210} \pm {\bf 2,340}</math></td></dl<>	$85,970 \pm 1,500$	$815.1 \pm 56.7$	${\bf 43,210} \pm {\bf 2,340}$
၀၀	$0.009075 \pm 0.002813$	<dl< td=""><td><dl< td=""><td><math>0.1777 \pm 0.0081</math></td><td><math>1.015 \pm 0.069</math></td></dl<></td></dl<>	<dl< td=""><td><math>0.1777 \pm 0.0081</math></td><td><math>1.015 \pm 0.069</math></td></dl<>	$0.1777 \pm 0.0081$	$1.015 \pm 0.069$
ö	$0.3442 \pm 0.0370$	<dl< td=""><td><dl< td=""><td><math>5.971 \pm 0.332</math></td><td><math display="block">\textbf{97.07} \pm \textbf{2.34}</math></td></dl<></td></dl<>	<dl< td=""><td><math>5.971 \pm 0.332</math></td><td><math display="block">\textbf{97.07} \pm \textbf{2.34}</math></td></dl<>	$5.971 \pm 0.332$	$\textbf{97.07} \pm \textbf{2.34}$
Cu	$3.853 \pm 1.519$	$0.7182 \pm 0.2176$	<dl< td=""><td><math>0.5375 \pm 0.0933</math></td><td><math display="block">\textbf{3,830} \pm \textbf{193}</math></td></dl<>	$0.5375 \pm 0.0933$	$\textbf{3,830} \pm \textbf{193}$
Fe	$26.21 \pm 0.83$	<dl< td=""><td><math>0.7639 \pm 0.2146</math></td><td><math>54.64 \pm 1.89</math></td><td><math>40.27 \pm 7.52</math></td></dl<>	$0.7639 \pm 0.2146$	$54.64 \pm 1.89$	$40.27 \pm 7.52$
Mg	$318.8 \pm 17.2$	$14.46 \pm 5.84$	$855.0 \pm 14.8$	$ 634,600 \pm 6,248 $	73,190 $\pm$ 3,860
Mn	$1.080 \pm 0.130$	$0.5131 \pm 0.2153$	$4.855 \pm 1.137$	$12.20 \pm 0.81$	$3,226 \pm 280$
Mo	$0.1232 \pm 0.0106$	<dl< td=""><td><dl< td=""><td><math>2.154 \pm 0.052</math></td><td><math>0.5949 \pm 0.0483</math></td></dl<></td></dl<>	<dl< td=""><td><math>2.154 \pm 0.052</math></td><td><math>0.5949 \pm 0.0483</math></td></dl<>	$2.154 \pm 0.052$	$0.5949 \pm 0.0483$
Ē	<dl< td=""><td><dl< td=""><td><dl< td=""><td><math>6.153 \pm 0.214</math></td><td><math>0.4588 \pm 0.1118</math></td></dl<></td></dl<></td></dl<>	<dl< td=""><td><dl< td=""><td><math>6.153 \pm 0.214</math></td><td><math>0.4588 \pm 0.1118</math></td></dl<></td></dl<>	<dl< td=""><td><math>6.153 \pm 0.214</math></td><td><math>0.4588 \pm 0.1118</math></td></dl<>	$6.153 \pm 0.214$	$0.4588 \pm 0.1118$
Se	<dl< td=""><td><math>0.01464 \pm 0.00514</math></td><td><dl< td=""><td><math>0.05102 \pm 0.00408</math></td><td><math display="block">\textbf{51.20} \pm \textbf{3.05}</math></td></dl<></td></dl<>	$0.01464 \pm 0.00514$	<dl< td=""><td><math>0.05102 \pm 0.00408</math></td><td><math display="block">\textbf{51.20} \pm \textbf{3.05}</math></td></dl<>	$0.05102 \pm 0.00408$	$\textbf{51.20} \pm \textbf{3.05}$
Zn	79,120 $\pm$ 1,710	$25.800 \pm 13.580$	<dl< td=""><td><dl< td=""><td><b>29,070</b> <math>\pm</math> <b>1,380</b></td></dl<></td></dl<>	<dl< td=""><td><b>29,070</b> <math>\pm</math> <b>1,380</b></td></dl<>	<b>29,070</b> $\pm$ <b>1,380</b>

<DL: Below detection limit. Elements in bold are not considered elemental contamination due to their inclusion in the ingredient list.



Concentrations ( $\mu$ g g<sup>-1</sup>) of Xenobiotic and Essential Elements Determined in Dietary Supplements from Company 3 (N = 12,95% Cls) TABLE 8.

	Zinc Supplement	Essential Fatty Acids Supplement	Calcium Supplement	Magnesium Supplement	Magnesium Supplement Multivitamin Supplement	Probiotic Supplement
Xenobio	Xenobiotic Elements					
₹	$15.44 \pm 6.29$	<dl< td=""><td><math>26.09 \pm 14.17</math></td><td><math display="block">31.93\pm4.82</math></td><td><math display="block">20.88\pm2.89</math></td><td><dl< td=""></dl<></td></dl<>	$26.09 \pm 14.17$	$31.93\pm4.82$	$20.88\pm2.89$	<dl< td=""></dl<>
As	$0.01153 \pm 0.00367$	<dl< td=""><td><math>0.02067 \pm 0.00647</math></td><td><math>0.4323 \pm 0.0221</math></td><td><math>0.5132 \pm 0.0238</math></td><td><math>0.04750 \pm 0.00780</math></td></dl<>	$0.02067 \pm 0.00647$	$0.4323 \pm 0.0221$	$0.5132 \pm 0.0238$	$0.04750 \pm 0.00780$
8	$0.03946 \pm 0.00618$	<dl< td=""><td><math>0.02777 \pm 0.00334</math></td><td><math>0.06962 \pm 0.00489</math></td><td><math>0.05674 \pm 0.00992</math></td><td><dl< td=""></dl<></td></dl<>	$0.02777 \pm 0.00334$	$0.06962 \pm 0.00489$	$0.05674 \pm 0.00992$	<dl< td=""></dl<>
βĤ	<dl< td=""><td><dl< td=""><td><dl< td=""><td><dl< td=""><td><dl< td=""><td><dl< td=""></dl<></td></dl<></td></dl<></td></dl<></td></dl<></td></dl<>	<dl< td=""><td><dl< td=""><td><dl< td=""><td><dl< td=""><td><dl< td=""></dl<></td></dl<></td></dl<></td></dl<></td></dl<>	<dl< td=""><td><dl< td=""><td><dl< td=""><td><dl< td=""></dl<></td></dl<></td></dl<></td></dl<>	<dl< td=""><td><dl< td=""><td><dl< td=""></dl<></td></dl<></td></dl<>	<dl< td=""><td><dl< td=""></dl<></td></dl<>	<dl< td=""></dl<>
Po	$0.4092 \pm 0.0074$	<dl< td=""><td><math>0.09268 \pm 0.00328</math></td><td><math>0.02407 \pm 0.00081</math></td><td><math>0.04092 \pm 0.00155</math></td><td><math>0.05856 \pm 0.00201</math></td></dl<>	$0.09268 \pm 0.00328$	$0.02407 \pm 0.00081$	$0.04092 \pm 0.00155$	$0.05856 \pm 0.00201$
Sp	<dl< td=""><td><dl< td=""><td><dl< td=""><td><math>0.02754 \pm 0.00169</math></td><td><dl< td=""><td><dl< td=""></dl<></td></dl<></td></dl<></td></dl<></td></dl<>	<dl< td=""><td><dl< td=""><td><math>0.02754 \pm 0.00169</math></td><td><dl< td=""><td><dl< td=""></dl<></td></dl<></td></dl<></td></dl<>	<dl< td=""><td><math>0.02754 \pm 0.00169</math></td><td><dl< td=""><td><dl< td=""></dl<></td></dl<></td></dl<>	$0.02754 \pm 0.00169$	<dl< td=""><td><dl< td=""></dl<></td></dl<>	<dl< td=""></dl<>
Sn	$0.1041 \pm 0.0131$	<dl< td=""><td><dl< td=""><td><math>0.01445 \pm 0.00362</math></td><td><math>0.4440 \pm 0.0260</math></td><td><math>0.02863 \pm 0.01034</math></td></dl<></td></dl<>	<dl< td=""><td><math>0.01445 \pm 0.00362</math></td><td><math>0.4440 \pm 0.0260</math></td><td><math>0.02863 \pm 0.01034</math></td></dl<>	$0.01445 \pm 0.00362$	$0.4440 \pm 0.0260$	$0.02863 \pm 0.01034$
⊃	$0.0007583 \pm 0.0001011$	< <b>D</b> L	$0.08226 \pm 0.00141$	$0.5130 \pm 0.0052$	$0.1032 \pm 0.0017$	$0.02199 \pm 0.00032$
>	$0.01963 \pm 0.00728$	<dl< td=""><td><math>0.3309 \pm 0.0123</math></td><td><math>2.833 \pm 0.062</math></td><td><math>11.66 \pm 3.27</math></td><td><math>0.01567 \pm 0.00683</math></td></dl<>	$0.3309 \pm 0.0123$	$2.833 \pm 0.062$	$11.66 \pm 3.27$	$0.01567 \pm 0.00683$
Essentia	Essential Elements					
S	$87.87 \pm 2.78$	<dl< td=""><td><math display="block">{\bf 299,000 \pm 79,600}</math></td><td>1, <math>150 \pm 80</math></td><td><math display="block">\textbf{43,000} \pm \textbf{740}</math></td><td><math>614.9 \pm 235.3</math></td></dl<>	${\bf 299,000 \pm 79,600}$	1, $150 \pm 80$	$\textbf{43,000} \pm \textbf{740}$	$614.9 \pm 235.3$
රි	<dl< td=""><td>&lt; DL</td><td><dl< td=""><td><math>0.02499 \pm 0.00344</math></td><td><math>0.4791 \pm 0.0120</math></td><td><math>0.07519 \pm 0.00497</math></td></dl<></td></dl<>	< DL	<dl< td=""><td><math>0.02499 \pm 0.00344</math></td><td><math>0.4791 \pm 0.0120</math></td><td><math>0.07519 \pm 0.00497</math></td></dl<>	$0.02499 \pm 0.00344$	$0.4791 \pm 0.0120$	$0.07519 \pm 0.00497$
ဝဲ	<dl< td=""><td>&lt;<b>D</b>L</td><td><dl< td=""><td><math>1.530 \pm 0.091</math></td><td><math display="block">\textbf{27.21} \pm \textbf{0.68}</math></td><td>&lt;<b>D</b>F</td></dl<></td></dl<>	< <b>D</b> L	<dl< td=""><td><math>1.530 \pm 0.091</math></td><td><math display="block">\textbf{27.21} \pm \textbf{0.68}</math></td><td>&lt;<b>D</b>F</td></dl<>	$1.530 \pm 0.091$	$\textbf{27.21} \pm \textbf{0.68}$	< <b>D</b> F
3	$0.1921 \pm 0.0192$	<dl< td=""><td><math>0.1949 \pm 0.0535</math></td><td><math>0.5806 \pm 0.0259</math></td><td><math>0.5046 \pm 0.0392</math></td><td><math>0.9184 \pm 0.0264</math></td></dl<>	$0.1949 \pm 0.0535$	$0.5806 \pm 0.0259$	$0.5046 \pm 0.0392$	$0.9184 \pm 0.0264$
æ	< <b>D</b> L	< DL	$21.18 \pm 3.19$	$44.97 \pm 1.02$	$43.61 \pm 2.28$	$3.924 \pm 0.802$
Mg	$33.21 \pm 3.68$	$2.480 \pm 0.326$	2, 423 $\pm$ 41	179,400 $\pm$ 3,200	$63,310 \pm 1,830$	$708.5 \pm 16.6$
Mn	$1.206 \pm 0.059$	< DL	$3.539 \pm 0.090$	$7.663 \pm 0.246$	$\textbf{360.8} \pm \textbf{4.8}$	$293.5\pm4.4$
Mo	<dl< td=""><td>&lt; DL</td><td><math>0.09942 \pm 0.00604</math></td><td><math>1.448 \pm 0.065</math></td><td><math display="block">\textbf{18.30} \pm \textbf{0.66}</math></td><td><math>0.08321 \pm 0.02328</math></td></dl<>	< DL	$0.09942 \pm 0.00604$	$1.448 \pm 0.065$	$\textbf{18.30} \pm \textbf{0.66}$	$0.08321 \pm 0.02328$
Z	<dl< td=""><td><dl< td=""><td><dl< td=""><td><math>0.8782 \pm 0.0715</math></td><td><math>0.4698 \pm 0.0798</math></td><td>&lt;<b>D</b>L</td></dl<></td></dl<></td></dl<>	<dl< td=""><td><dl< td=""><td><math>0.8782 \pm 0.0715</math></td><td><math>0.4698 \pm 0.0798</math></td><td>&lt;<b>D</b>L</td></dl<></td></dl<>	<dl< td=""><td><math>0.8782 \pm 0.0715</math></td><td><math>0.4698 \pm 0.0798</math></td><td>&lt;<b>D</b>L</td></dl<>	$0.8782 \pm 0.0715$	$0.4698 \pm 0.0798$	< <b>D</b> L
Se	<dl< td=""><td><math>0.008819 \pm 0.001073</math></td><td><math>0.02959 \pm 0.00397</math></td><td><dl< td=""><td><math display="block">\textbf{50.80} \pm \textbf{2.73}</math></td><td><math>0.01847 \pm 0.00143</math></td></dl<></td></dl<>	$0.008819 \pm 0.001073$	$0.02959 \pm 0.00397$	<dl< td=""><td><math display="block">\textbf{50.80} \pm \textbf{2.73}</math></td><td><math>0.01847 \pm 0.00143</math></td></dl<>	$\textbf{50.80} \pm \textbf{2.73}$	$0.01847 \pm 0.00143$
Zu	$\textbf{72.08} \pm \textbf{1.02}$	<dl< td=""><td><math>38.69 \pm 1.67</math></td><td><math>4.337 \pm 0.233</math></td><td><math display="block">\textbf{2,839} \pm \textbf{123}</math></td><td><math>24.62 \pm 2.53</math></td></dl<>	$38.69 \pm 1.67$	$4.337 \pm 0.233$	$\textbf{2,839} \pm \textbf{123}$	$24.62 \pm 2.53$

<DL: Below detection limit. Elements in bold are not considered elemental contamination due to their inclusion in the ingredient list.



Concentrations ( $\mu$ g g<sup>-1</sup>) of Xenobiotic and Essential Elements Determined in Dietary Supplements from Company 4 (N = 12, 95% Cls) TABLE 9.

Zinc Supplement		Supplement	Calcium Supplement	Magnesium Supplement	Magnesium Supplement Multivitamin Supplement	Probiotic Supplement
Xenobiotic Elements	8.8	3.424±1.208	$91.74 \pm 3.42$	<dl< td=""><td>31.22 ± 3.050</td><td><dl< td=""></dl<></td></dl<>	31.22 ± 3.050	<dl< td=""></dl<>
^ □^	0.0075	$0.007569 \pm 0.002793$	$0.2096 \pm 0.0125$	$0.01163 \pm 0.00413$	$0.09037 \pm 0.00767$	$0.02217 \pm 0.00572$
<dl></dl>		<dl< td=""><td><math>0.1194 \pm 0.0077</math></td><td><dl< td=""><td><math>0.03310 \pm 0.00605</math></td><td><dl< td=""></dl<></td></dl<></td></dl<>	$0.1194 \pm 0.0077$	<dl< td=""><td><math>0.03310 \pm 0.00605</math></td><td><dl< td=""></dl<></td></dl<>	$0.03310 \pm 0.00605$	<dl< td=""></dl<>
<sup>∨</sup> DL		<dl <<="" td=""><td><dl< td=""><td><dl< td=""><td><dl< td=""><td>^DL</td></dl<></td></dl<></td></dl<></td></dl>	<dl< td=""><td><dl< td=""><td><dl< td=""><td>^DL</td></dl<></td></dl<></td></dl<>	<dl< td=""><td><dl< td=""><td>^DL</td></dl<></td></dl<>	<dl< td=""><td>^DL</td></dl<>	^DL
$0.003119 \pm 0.000568$	000568	<dl< td=""><td><math>0.05814 \pm 0.00626</math></td><td><dl< td=""><td><math>0.06468 \pm 0.00132</math></td><td><dl< td=""></dl<></td></dl<></td></dl<>	$0.05814 \pm 0.00626$	<dl< td=""><td><math>0.06468 \pm 0.00132</math></td><td><dl< td=""></dl<></td></dl<>	$0.06468 \pm 0.00132$	<dl< td=""></dl<>
$14.85 \pm 0.15$	15	<dl< td=""><td>OL</td><td><dl< td=""><td><dl< td=""><td><dl< td=""></dl<></td></dl<></td></dl<></td></dl<>	OL	<dl< td=""><td><dl< td=""><td><dl< td=""></dl<></td></dl<></td></dl<>	<dl< td=""><td><dl< td=""></dl<></td></dl<>	<dl< td=""></dl<>
^ PL		<dl <<="" td=""><td><dl< td=""><td><math>0.03139 \pm 0.00902</math></td><td><dl< td=""><td>^DL</td></dl<></td></dl<></td></dl>	<dl< td=""><td><math>0.03139 \pm 0.00902</math></td><td><dl< td=""><td>^DL</td></dl<></td></dl<>	$0.03139 \pm 0.00902$	<dl< td=""><td>^DL</td></dl<>	^DL
^ PL		<dl< td=""><td><math>0.3197 \pm 0.0057</math></td><td><math>0.004984 \pm 0.000197</math></td><td><math>0.08917 \pm 0.00175</math></td><td><math>0.01587 \pm 0.00036</math></td></dl<>	$0.3197 \pm 0.0057$	$0.004984 \pm 0.000197$	$0.08917 \pm 0.00175$	$0.01587 \pm 0.00036$
<dl< td=""><td>0.0072</td><td><math>0.007211 \pm 0.002768</math></td><td><math>1.719 \pm 0.035</math></td><td><math>0.8498 \pm 0.0284</math></td><td><math>0.5660 \pm 0.0126</math></td><td><math display="block">0.01093 \pm 0.00424</math></td></dl<>	0.0072	$0.007211 \pm 0.002768$	$1.719 \pm 0.035$	$0.8498 \pm 0.0284$	$0.5660 \pm 0.0126$	$0.01093 \pm 0.00424$
Essential Elements						
<sup>∨</sup> DL	6.1	$153 \pm 1.613$	${\bf 242,900} \pm {\bf 4,100}$	$198.7 \pm 24.2$	$  16,760 \pm 1,230 $	$73.11 \pm 18.88$
^DF		<dl< td=""><td><math>0.03813 \pm 0.00433</math></td><td><math>0.05813 \pm 0.00423</math></td><td><math>0.2548 \pm 0.0055</math></td><td><math>0.05438 \pm 0.00386</math></td></dl<>	$0.03813 \pm 0.00433$	$0.05813 \pm 0.00423$	$0.2548 \pm 0.0055$	$0.05438 \pm 0.00386$
^DL		<dl< td=""><td><math>2.430 \pm 0.201</math></td><td><math>2.813 \pm 0.096</math></td><td><math display="block">\textbf{14.32} \pm \textbf{0.42}</math></td><td><dl< td=""></dl<></td></dl<>	$2.430 \pm 0.201$	$2.813 \pm 0.096$	$\textbf{14.32} \pm \textbf{0.42}$	<dl< td=""></dl<>
^DL		<dl< td=""><td><math>0.2673 \pm 0.0250</math></td><td><math>0.1516 \pm 0.0153</math></td><td><math>0.08997 \pm 0.01458</math></td><td><math>0.3461 \pm 0.0134</math></td></dl<>	$0.2673 \pm 0.0250$	$0.1516 \pm 0.0153$	$0.08997 \pm 0.01458$	$0.3461 \pm 0.0134$
< <u>0</u> > DL		<dl< td=""><td><math>133.2 \pm 2.8</math></td><td><math>29.04 \pm 1.26</math></td><td><math>39.62 \pm 1.22</math></td><td><dl< td=""></dl<></td></dl<>	$133.2 \pm 2.8$	$29.04 \pm 1.26$	$39.62 \pm 1.22$	<dl< td=""></dl<>
< <b>D</b>	5.5	$118 \pm 0.709$	1, $459 \pm 23$	196,500 $\pm$ 4,300	$\textbf{19,020} \pm \textbf{970}$	1, 440 $\pm$ 30
< <b>D</b>		<dl< td=""><td><math>44.68 \pm 1.31</math></td><td><math>4.370 \pm 0.147</math></td><td><math display="block">\textbf{920.2} \pm \textbf{60.5}</math></td><td><math>221.0 \pm 3.5</math></td></dl<>	$44.68 \pm 1.31$	$4.370 \pm 0.147$	$\textbf{920.2} \pm \textbf{60.5}$	$221.0 \pm 3.5$
<sup>∨</sup> DL		<dl< td=""><td>&lt; DL</td><td><math>0.6047 \pm 0.0329</math></td><td><math>9.774 \pm 0.742</math></td><td><dl< td=""></dl<></td></dl<>	< DL	$0.6047 \pm 0.0329$	$9.774 \pm 0.742$	<dl< td=""></dl<>
^DL		<dl< td=""><td>&lt; DL</td><td><math>1.164 \pm 0.062</math></td><td><dl< td=""><td>~DL</td></dl<></td></dl<>	< DL	$1.164 \pm 0.062$	<dl< td=""><td>~DL</td></dl<>	~DL
<dl></dl>	0.024	$0.02422 \pm 0.00209$	$0.006817 \pm 0.001101$	$0.01070 \pm 0.00145$	$\textbf{12.44} \pm \textbf{1.39}$	$0.01334 \pm 0.00159$
$\textbf{2,829} \pm \textbf{81}$	31	<dl< td=""><td><math>4.298 \pm 0.260</math></td><td><dl< td=""><td><b>6,323</b> <math>\pm</math> 220</td><td><math>10.65 \pm 0.23</math></td></dl<></td></dl<>	$4.298 \pm 0.260$	<dl< td=""><td><b>6,323</b> <math>\pm</math> 220</td><td><math>10.65 \pm 0.23</math></td></dl<>	<b>6,323</b> $\pm$ 220	$10.65 \pm 0.23$

<DL: Below detection limit. Elements in bold are not considered elemental contamination due to their inclusion in the ingredient list.



tion of each more than 15  $\mu$ g g<sup>-1</sup>. Except for the essential fatty acids supplement, Pb was quantified in all of the samples with the zinc supplement having the highest concentration at nearly 0.41  $\mu$ g g<sup>-1</sup>. Of all the supplements tested, the multivitamin supplement had the highest concentration of As at 0.513  $\mu$ g g<sup>-1</sup>, and the magnesium supplement contained nearly as much As at 0.432  $\mu$ g g<sup>-1</sup> along with  $0.444 \ \mu g \ g^{-1}$  of Sn. In both calcium and zinc supplements, elemental contamination from V, As, and Cd was found. The probiotic had low concentrations of V, As, Sn, and U. The essential fatty acids supplement was clear of contamination from xenobiotic elements with values below the respective detection limits.

Similar to both Companies 2 and 3, Al was a common contaminant among products of Company 4, as it was found at nearly 100  $\mu$ g g<sup>-1</sup> in the calcium supplement, over 30  $\mu g$  g<sup>-1</sup> in the multivitamin supplement, and above 3  $\mu g$  g<sup>-1</sup> in the essential fatty acids supplement (Table 9). The zinc supplement contained nearly 15  $\mu$ g g<sup>-1</sup> of Sb. In the calcium supplement, V was quantified at over 1.7  $\mu$ g g<sup>-1</sup> along with determining the presence of As, Cd, and U. Vanadium was found at 0.85  $\mu$ g g<sup>-1</sup> in the magnesium supplement in addition to As, Sn, and U. An elevated V amount was also observed in the multivitamin supplement at  $0.566 \ \mu g \ g^{-1}$  along with As, Cd, Pb, and U. As and U were found in the probiotic supplement.

### Essential Elements

In addition to the presence of xenobiotic elements, each dietary supplement from Company 1 had notable contamination by essential elements that were not listed as ingredients on the dietary supplement labels (Table 6). Overall, the calcium and magnesium supplements contained higher quantities of unlabeled essential elements compared with xenobiotic elements. In the calcium supplement, Mg was found at a concentration of over 5,000  $\mu$ g g<sup>-1</sup>, Zn at more than  $80~\mu$ g g<sup>-1</sup>, and Fe at  $8 \mu g g^{-1}$ . The amount of contamination in the magnesium supplement consisted of Ca quantified at 200  $\mu$ g g<sup>-1</sup> and Fe at about 15  $\mu$ g g<sup>-1</sup> along with the presence of Cr, manganese (Mn), and nickel (Ni). In contrast, the zinc supplement contained no detectable quantities of any essential element.

Similar to Company 1, essential elements not listed as ingredients were found at a high frequency in the dietary supplements from Company 2 (Table 7). Magnesium was quantified at 855  $\mu$ g g<sup>-1</sup> in the calcium supplement and 318  $\mu$ g g<sup>-1</sup> in the zinc supplement. Calcium was found at 815 and 547  $\mu$ g g<sup>-1</sup> in the magnesium and zinc supplements, respectively. Except for the essential fatty acids supplement, Fe was observed in all the samples with the highest concentration of 54  $\mu$ g g<sup>-1</sup> in the magnesium supplement. The magnesium supplement was also determined to contain Cr, Mn, molybdenum (Mo), and Ni at concentrations between 2.1 and  $12.2 \mu g g^{-1}$ . The calcium supplement had a Mn concentration of nearly 5  $\mu$ g g<sup>-1</sup>. Cobalt (Co), Mo, and Ni were detected in the multivitamin supplement. The zinc supplement frequently contained essential elements, notably Cu, Mn, and Cr. The presence of Zn and Mg was discovered in the essential fatty acids supplement.

The products from Company 3 also had contamination from unlabeled essential elements (Table 8). There was nearly  $300 \,\mu\mathrm{g}\,\mathrm{g}^{-1}$  of Mn in the probiotic supplement. Magnesium was quantified at 2,500, 700, 33, and 2.5  $\mu$ g g<sup>-1</sup> in the calcium, probiotic, zinc, and the essential fatty acids supplements, respectively. In the magnesium and



zinc supplements, 1,150 and 88  $\mu$ g g<sup>-1</sup> Ca was measured, respectively. Zinc was also found in the calcium, magnesium, and probiotic supplements. Iron was detected in the magnesium, multivitamin, calcium, and probiotic supplements. Manganese was observed in the magnesium, calcium, and zinc supplements. Copper was found in the calcium, magnesium, zinc, and probiotic supplements. The magnesium supplement also had Cr, Ni, and Co. In addition, the multivitamin supplement contained Co and Ni.

As with the other three companies, the dietary supplements from Company 4 also contained contamination from essential elements (Table 9). Similar to Company 3, the probiotic supplement from Company 4 had a Mn concentration close to 250  $\mu$ g g<sup>-1</sup>. Magnesium was quantified at greater than 1,000  $\mu$ g g<sup>-1</sup> in the calcium and probiotic supplements. Calcium was measured at 198 and 73  $\mu$ g g<sup>-1</sup> in the magnesium and probiotic supplements. Of all the essential fatty acid supplements, the product from Company 4 was the most contaminated with essential elements, including Ca and Mg. In the probiotic and calcium supplements, Zn was quantified at 10 and 4  $\mu$ g g<sup>-1</sup>, respectively. Iron was found in the calcium, multivitamin, and magnesium supplements. In the calcium supplement, the Mn concentration was determined to be more than 40  $\mu$ g g<sup>-1</sup> with Cr at 2.4  $\mu$ g g<sup>-1</sup>. The magnesium supplement had Cr, Mn, and Ni at concentrations higher than 1  $\mu$ g g<sup>-1</sup> along with Mo. Except for the zinc and essential fatty acid supplements, all of the supplements from Company 4 contained Cu.

### DISCUSSION

The high fidelity results obtained in this study were accurate and precise based on validation with SRM 3280. The microwave-enhanced sample preparation provided effective, consistent, and efficient sample decomposition that combined with the robust ICP-MS analysis demonstrated the feasibility of using EPA methods to routinely determine both xenobiotic and essential elemental contamination in dietary supplements. Only two other studies have used SRM 3280 to validate the method protocols employed for multi-element quantification in dietary supplements or botanicals (Avula et al., 2010, 2011). This study demonstrated superior analytical and statistical procedures compared with the above-mentioned two previous studies by Avula et al. through increased quality control and inclusion of CIs for all sample results. Another important improvement in the current study was related to the preparation of the calibration standards. The previous studies used a calibration blank solution of 3.2% HNO<sub>3</sub> and 0.8% HCl (Avula et al., 2010, 2011). Although these studies produced accurate results for SRM 3280, this method is not optimal in providing uniform acid content in all of the digested samples and calibration solutions. However, the addition of appropriate amount of analytical blank solution to each calibration standard during the subsequent dilution as performed in the current study accomplished this goal.

This study confirmed the presence of contamination from either xenobiotic or essential elements in all the samples from all four suppliers. The xenobiotic elements quantified in this study included Pb, Cd, Sb, Sn, and As, with their concentrations ranging from below  $1 \mu g g^{-1}$  to  $\mu g g^{-1}$  levels. Only two supplements did not have xenobiotic elements detected in them, the zinc supplement from Company 1



and the essential fatty acid supplement from Company 2. Particularly alarming was the frequent presence of Al in  $\mu g g^{-1}$  quantities in half of the supplements tested. The most common sources of xenobiotic elemental contamination were V and U. Nearly two-thirds of the samples contained quantifiable amounts of Pb. Antimony was found in only one supplement at a concentration more than than 0.05  $\mu$ g g<sup>-1</sup>. One positive note was that the Hg concentrations were below the detection limit in all of the samples tested.

In addition to the presence of xenobiotic elements, most of the supplements contained essential elements such as Zn, Mg, Ca, Fe, Cu, Mo, and Se that were not listed as ingredients on the labels. This study identified a higher frequency of adulteration from essential elements compared with xenobiotic elements, with Ca being the most frequent contaminant. The zinc supplements from Companies 1 and 4 were the only two that did not contain quantifiable amounts of any extra essential element. Iron was almost always found in  $\mu g g^{-1}$  concentration as an impurity. The presence of either xenobiotic or essential elements in all of the dietary supplements raises the possibility of a systemic and persistent problem across the industry.

The contamination observed in this study was examined through a comparison of the four manufacturers. Company 3 had the highest occurrence of xenobiotic elemental content. Except for the essential fatty acid, each supplement contained at least five xenobiotic elements. It was also the only company that had four supplements with  $\mu g g^{-1}$  quantities of Al along with the two highest As concentrations. Their multivitamin supplement had the highest amount of As and V quantified in the study, while their zinc supplement had the second highest concentration of Pb. The elemental adulteration for Company 4 was nearly as frequent as Company 3. Their calcium supplement had the highest concentration of Al with the highest Sb concentration found in their zinc supplement. Company 1 had the least amount of xenobiotic content. However, its calcium supplement contained the highest concentration of Sn, and their zinc supplement had the second highest concentration of Cd. Both calcium and magnesium supplements from Company 1 had a similar pattern of xenobiotic contamination from Pb, Sn, U, and V. The company's use of the same ingredients in different amounts in their supplements may provide one explanation (Table 2). These supplement types from Companies 2 and 3 shared similar elemental adulteration but this was not observed in Company 4. Three of the five supplements from Company 2 had at least 10  $\mu$ g g<sup>-1</sup> of Al. The zinc supplement from Company 2 contained the highest concentration of Pb.

Some patterns appeared when the supplement types were compared. The essential fatty acids were the cleanest for xenobiotic elements as the products from Company 4 only contained over 3  $\mu$ g g<sup>-1</sup> of Al. The essential fatty acids were more likely to contain essential elements with Mg detected in all three samples analyzed. The zinc supplements were relatively clear of unlabeled essential elements, except for the one from Company 2 that contained Ca and Mg along with six other essential elements. The contamination frequency was less for the zinc supplement from Company 3, but it still included  $\mu g g^{-1}$  quantities of Ca and Mg. All of the calcium supplements were adulterated with at least 855  $\mu$ g g<sup>-1</sup> of Mg. The four magnesium supplements contained at least 198  $\mu$ g g<sup>-1</sup> of Ca and 15.8  $\mu$ g g<sup>-1</sup> of



Fe. Even though they had the most elements listed as ingredients, the multivitamins contained contamination from nearly all of the other essential elements, notably Fe at approximately 40  $\mu$ g g<sup>-1</sup> in all three. Both of the probiotics had over 200  $\mu$ g g<sup>-1</sup> of Mn. Certain bacteria are known to produce Mn, especially the Lactobacillus genus (Bomba et al., 2002), a common species in this type of supplement. The bacteria also need the other essential elements that were quantified in the probiotics for survival, especially Ca and Mg that were found in high concentrations. The media where the bacteria were grown is a possible contamination source.

The findings of this study are significant because they reveal the quality of dietary supplements in the current marketplace and the possible threat they pose to public health as their use increases globally. The long-term health effects from regular intake and exposure to elemental contaminants from dietary supplements are rarely studied and insufficiently understood. One of the less publicized aspects of these quality deficiencies is the presence of essential elements not mentioned on content labels, potentially negatively affecting supplement regimens designed to improve human health through nutrient support. With people taking multiple dietary supplements, potential negative synergistic effects are a realistic possibility and cumulative exposure to these xenobiotic elements may have long-term adverse health outcomes. In addition to the consequences from the cumulative exposure to these potentially harmful elements, over-supplementation can result from the presence of essential elements that were not listed on the label, as demonstrated by this study.

This study demonstrated the application of standardized EPA methods for the evaluation of elemental contamination and content in dietary supplements. The use of these methods resulted in finding xenobiotic and/or essential elemental contamination in all the dietary supplements analyzed. Their implementation in both sample preparation and analysis to determine elemental concentrations could serve as the foundation of a reliable, consistent, and efficient quality assurance program for the dietary supplement industry and regulatory agencies. The overwhelming trend that emerged from this study was the widespread contamination in the products from the different suppliers. The elemental adulteration observed could result from the manufacturing process or the raw ingredients that comprise the supplements. It is difficult to completely know their sources, but many ingredients, such as Stevia, come from countries where regulation is not as stringent as in the United States. It is clear from this study that in order to provide the highest quality of dietary supplements possible, all final products and ingredients must be analyzed for elemental contamination to ultimately meet the intent of the GMP standards of the Dietary Supplement Health and Education Act of 1994. Further studies are necessary to evaluate the type, level, and frequency of dietary supplement contamination throughout the industry.

### ACKNOWLEDGMENTS

The authors would like to acknowledge Dr. Gerra Bosco for her assistance in the preparation of this manuscript. Dan Bodnar, Dave Hardesty, Lance Crosby, and Ben Lauterbach were responsible for instrument maintenance. Duquesne



University and Applied Isotope Technologies (AIT) supplied the resources to perform the study. The authors thank Agilent Technologies and Milestone for instrumentation support.

### ABOUT THE AUTHORS

**Dr. Zinn** received his Ph.D. in Analytical Chemistry from Duquesne University in 2014. His dissertation focused on the evaluation and assessment of autism spectrum disorders through the application of mass spectrometric measurement techniques. **Dr. Rahman** is a senior chemist at BASF. He previously worked as the R&D Manager at AIT where he developed and performed both elemental and speciation methods using ICP-MS. **Dr. Faber** is a developmental pediatrician who works at The Children.s Institute of Pittsburgh. He provides outpatient services for children with severe developmental and behavioral disorders along with conducting clinical research involving children with autism spectrum disorders. **Dr. Wolle** is a Research Fellow at the US Food and Drug Administration. He previously worked as Post-Doctoral Fellow at Duquesne University during which he worked to develop speciation methods using ICP-MS. Mr. Pamuku is the CEO of AIT, a company that provides elemental, speciation, and molecular measurement solutions using advanced measurement techniques. Dr. Kingston is a Professor of Analytical Chemistry at Duquesne University in the Department of Chemistry and Biochemistry. His research focuses on the development and application of mass spectrometric methods in environmental human health, most notably Speciated Isotope Dilution Mass Spectrometry.

### REFERENCES

- Aldosary BM, Sutter ME, Schwartz M, Morgan BW. Case series of selenium toxicity from a nutritional supplement. Clin Toxicol. 2012;50(1):57–64.
- Ashar BH, Rice TN, Sisson SD. Physicians' understanding of the regulation of dietary supplements. Arch Intern Med. 2007;167(9):966-969.
- Avula B, Wang YH, Smillie TJ, Duzgoren-Aydin NS, Khan IA. Quantitative determination of multiple elements in botanicals and dietary supplements using ICP-MS. J Agr Food Chem. 2010;58(16):8887-8894.
- Avula B, Wang YH, Duzgoren-Aydin NS, Khan IA. Inorganic elemental conpositions of commercial multivitamin/mineral dietary supplements: application of collision/reaction cell inductively coupled-mass spectroscopy. Food Chem. 2011;127(1):54–62.
- Bomba A, Nemcova R, Mudronova D, Guba P. The possibilities of potentiating the efficacy of probiotics. Trends Food Sci Technol. 2002;13(4):121–126.
- Bourgoin BP, Evans DR, Cornett JR, Lingard SM, Quattrone AJ. Lead content in 70 brands of dietary calcium supplements. Am J Public Health. 1993;83(3):1155–1160.
- Brewer GJ. The risks of copper toxicity contributing to cognitive decline in the aging population and to Alzheimer's disease. J Am Coll Nutr. 2009;28(3):238-242.
- Brewer GJ. Issues raised involving the copper hypothesis in the causation of Alzheimer's disease. Int J Alzheimer's Dis. 2011;537528.
- Chen WY, Wang YC, Kuo M-S. Determination of total mercury and methylmercury in human hair by graphite-furnace atomic absorption spectrophotometry using 2,3-Dimercaptopropane-1-sulfonate as a complexing agent. Anal Sci. 2002;18(3):255–260.
- Chuang IC, Chen KS, Huang Y-L, Lee P-N, Lin T-H. Determination of trace elements in some natural drugs by atomic absorption spectrometry. Biol Trace Elem Res. 2000;76(3):235-244.



- Cohen PA. American roulette-contaminated dietary supplements. The New Engl J Med. 2009;361(16):1523-1525.
- Crowley, R, FitzGerald LH. The impact of cGMP compliance on consumer confidence in dietary supplement products. Toxicology. 2006;221(1):9–16.
- Dolan SP, Nortrup DA, Bolger PM, Capar SG. Analysis of dietary supplements for arsenic, cadmium, mercury, and lead using inductively coupled plasma mass spectrometry. J Agr Food Chem. 2003;51(5):1307–1312.
- Driscoll MS, Kwon EKM, Skupsky H, Kwon S-Y, Grant-Kels JM. Nutrition and the deleterious side effects of nutritional supplements. Clin Dermatol. 2010;28(4):371–379.
- Fairweather-Tait SJ, Harvey LJ, Collings R. Risk-benefit analysis of mineral intakes: case studies on copper and iron. Proc Nutr Soc. 2011;70(1):1–9.
- Fan AM, Kizer KW. Selenium-nutritional, toxicological, and clinical aspects. The West J Med. 1990:153(2):160-167.
- Garcia-Rico L, Leyva-Perez J, Jara-Marini ME. Content and daily intake of copper, zinc, lead, cadmium, and mercury from dietary supplements in Mexico. Food Chem Toxicol. 2007;45(9):1599–1605.
- Geyer H, Parr MK, Koehler K, Mareck U, Schanzer W, Thevis M. Nutritional supplements crosscontaminated and faked with doping substances. J Mass Spectrom. 2008;43(7):892–902.
- Igic PG, Lee E, Harper W, Roach KW. Toxic effects associated with consumption of zinc. Mayo Clin Proc. 2002;77(7):713–716.
- Jomova K, Jenisova Z, Feszterova M, Baros S, Liska J, Hudecova D, Rhodes CJ, Valko M. Arsenic: toxicity, oxidative stress and human disease. J Appl Toxicol. 2011;31(2):95–107.
- Kawahara, M, Kato-Negishi M. Link between aluminum and the pathogenesis of Alzheimer's disease: the integration of the aluminum and amyloid cascade hypotheses. Int J Alzheimer's Dis. 2011:2011.
- Kingston HMS, Walter PJ, Chalk S, Lorentzen E, Link D. Environmental microwave sample preparation: fundamentals, methods, and applications. In: HMS Kingston and SJ Haswell (Eds.), Microwave-enhanced chemistry: fundamentals, sample preparation, and applications, pp 223–349, Washington DC: American Chemical Society, 1997.
- Linder MC, Hazegh-Azam M. Copper biochemistry and molecular biology. Am J Clin Nutr. 1996;63(5):797S-811S.
- Lopez de Romana D, Olivares M, Uauy R, Araya M. Risks and benefits of copper in light of new insights of copper homeostasis. J Trace Elem Med Biol. 2011;25(1):3–13.
- MacFarquhar JK, Broussard DL, Melstrom P, Hutchinson R, Wolkin A, Martin C, Burk RF, Dunn JR, Green AL, Hammond R, Schaffner W, Jones TF. Acute selenium toxicity associated with a dietary supplement. Arch Intern Med. 2010;170(3):256–261.
- Marrero J, Rebagliati RJ, Leiva E, Londonio A, Smichowski P. Inductively coupled plasma optical emission spectrometric determination of fifteen elements in dietary supplements: are the concentrations declared in the labels accurate? Microchem J. 2013;108:81–86.
- Maughan RJ. Contamination of dietary supplements and positive drug tests in sport. J Sports Sci. 2005;23(9):883–889.
- Maughan RJ, Greenhaff PL, Hespel P. Dietary supplements for athletes: emerging trends and recurring themes. J Sports Sci. 2011;29(S1):S57–S66.
- McCormick DB. Vitamin/mineral supplements: of questionable benefit for the general population. Nutr Rev. 2010;68(4):207-213.
- Mindak WR, Cheng J, Canas BJ, Bolger PM. Lead in women's and children's vitamins. J Agr Food Chem. 2008;56(16):6892–6896.
- Neustadt, J, Pieczenik S. Heavy-metal toxicity-with emphasis on mercury. Integr Med. 2007;6(2):26-32.
- Nuttall KL. Evaluating selenium poisoning. Ann Clin Lab Sci. 2006;36(4):409–420.
- Obi E, Akunyili DN, Ekpo B, Orisakwe OE. Heavy metal hazards of Nigerian herbal remedies. Sci Total Environ. 2006;369(1–3):35–41.
- Office of Environmental Health Hazard Assessment. State of California Environmental Protection Agency Office of Environmental Health Hazard Assessment Safe Drinking



Water and Toxic Enforcement Act of 1986. Retrieved December 18, 2013, from http://www.oehha.org/prop65.html.

- Porea TJ, Belmont JW, Mahoney Jr. DH. Zinc-induced anemia and neutropenia in an adolescent. The J Pediatr. 2000;136(5):688-690.
- Radimer K, Bindewald B, Hughes J, Ervin B, Swanson C, Piccano MF. Dietary supplement use by US adults: data from the National Health and Nutrition Examination Survey, 1999-2000. Am J Epidemiol. 2004;160(4):339–349.
- Raghunath R, Tripathi RM, Suseela B, Bhalke S, Shukla VK, Puranik VD. Dietary intake of metals by Mumbai adult population. Sci Total Environ. 2006;356(1-3):62-68.
- Raman P, Patino LC, Nair MG. Evaluation of metal and microbial contamination in botanical supplements. J Agr Food Chem. 2004;52(26):7822–7827.
- Saldanha LG, Betz JM, Coates PM. Development of the analytical methods and reference materials program for dietary supplement at the National Institute of Health. J AOAC Int. 2004;87(1):162-165.
- Salzman MB, Smith EM, Koo C. Excessive oral zinc supplementation. J Pediatr Hematol/Oncol. 2002;24(7):582–584.
- Sargent JD. The role of nutrition in the prevention of lead poisoning in children. Pediatr Ann 1994;23(11):636-642.
- Scelfo GM, Flegal AR. Lead in calcium supplements. Environ Health Perspect. 2000;108(4):309-313.
- Shankar AH, Prasad AS. Zinc and immune function: the biological basis of altered resistance to infection. Am J Clin Nutr. 1998;68(2 Suppl):447S-463S.
- Skoog DA, Holler FJ, Crouch SR. Principles of instrumental analysis. Australia: Thomson Brooks/Cole, 2007.
- Sutter ME, Thomas JD, Brown J, Morgan B. Selenium toxicity: a case of selenosis caused by a nutritional supplement. Ann Intern Med. 2008;148(12):970–971.
- Taylor CL. Regulatory frameworks for functional foods and dietary supplements. Nutrit Rev. 2004;62(2):55–59.
- Thomas R. Practical guide to ICP-MS a tutorial for beginners. Boca Raton, FL: CRC Press, 2008. Tumir H, Bosnir J, Vedrina-Dragojevic I, Dragun Z, Tomic S, Puntaric D, Jurak G. Monitoring of metal and metalloid content in dietary supplements on the Croatian market. Food Contr. 2010;21:885-889.
- US Environmental Protection Agency (EPA). US EPA method 6020B: inductively coupled plasma-mass spectrometry in test methods for evaluating solid waste, physical/chemical methods SW-846. Washington, DC: Government Printing Office, 2013.
- US Environmental Protection Agency (EPA). National Primary Drinking Water Regulations. Retrieved December 18, 2013, from http://water.epa.gov/drink/contaminants/index.cfm.
- United States Pharmacopeia. <232> Elemental impurities—limits. Second Supplement to USP 35-NF 30, United States Pharmacopeial Forum. 2012a;40(2):1-3
- United States Pharmacopeia. <2232> Elemental contaminants in dietary supplements. Pharm Forum. 2012b;38(3):1–3
- Woo JJY. Adverse event monitoring and multivitamin-multimineral dietary supplements. The Am J Clin Nutr. 2007;85(1):323S-324S.

