Final Report

Mammalian Erythrocyte Micronucleus Test of MWCNT in ICR mice

(Study No.: GT13-00019)

November 2014



Bioconvergence Technology Laboratory

Statement

Study No. : GT13-00019

Title: Mammalian Erythrocyte Micronucleus Test of MWCNT in ICR mice

This final report was written in Korean and translated into English.

This study has been performed in compliance with the principles of Good Laboratory Practices and test guidelines in following documents.

- National Institute of Environment Research (NIER) [Notice No. 2012-23, (revised 22 August 2012)].
- 2. OECD Guidelines for the Testing of Chemicals No. 474 "Mammalian Erythrocyte Micronucleus test" (Adopted: 21st July 1997).

The stated object in study protocol was achieved and there were no significant deviations from the aforementioned regulations that affected quality or integrity of the study. Therefore the justification of all data in this study was confirmed. The information of the test substance was written from the document that study sponsor provided.

7/6	14 November 2014
Jin-Sik Kim	Date

Study Director

Bioconvergence Technology Laboratory, Korea Conformity Laboratories

Jin-Kyu Lee Nov. (4, vol4)

Date

Management

Bioconvergence Technology Laboratory, Korea Conformity Laboratories

QUALITY ASSURANCE STATEMENT

Study No.: GT13-00019

Title: Mammalian Erythrocyte Micronucleus Test of MWCNT in ICR

mice

This study was subject to audit by the independent Quality Assurance Unit of KCL as indicated below. The findings of each audit were reported to the study director and management as prescribed by Standard Operating Procedures.

The final report audit was designed to confirm that as far as can be reasonably established the methods described and results incorporated in the final report accurately reflect the raw data produced during the study.

Audit phases and dates reported to the responsible personnel were as indicated below and these were based upon the audit records.

Phase Inspected	Date	Reports to Study Director	Reports to Management	
Study Plan	2013. 03. 18	2013. 03. 18	2013. 03. 18	
Animal receipt	2013, 05, 14	2013. 05. 14	2013. 05. 14	
Storage of test substance and vehicle	2013. 04. 01	2013. 04. 01	2013. 04. 01	
Preparation of test substance	2013, 05, 21	2013. 05. 21	2013. 05. 21	
Animal care and Administration	2013. 05. 21	2013; 05, 21	2013. 05, 21	
Clinical sign	2013. 05, 21	2013. 05. 21	2013. 05. 21	
Necropsy and Preparation of specimen	2013. 05. 22	2013. 05. 22	2013, 05, 22	
Observation	2013. 05. 22	2013. 05. 22	2013. 05. 22	
Raw data	2013. 08. 23	2013. 08. 23	2013. 08. 23	
Final report	2013, 08, 23	2013, 08, 23	2013. 08. 23	

QA director :

Song, Kyung Seuk Ph.D. Auditor, Quality Assurance Date 2013. 08. 23

* signed original

Study Personnel

Techincal Assistant Hyo-Jin Joo * Date 23 August 2013 23 August 2013 Sample preparation Jae-Hyuck Sung * **Date Animal Experiment** Hye-Jin Kim * **Date** 23 August 2013 Archiving of documents Hyo-Dong Kim * Date 23 August 2013

Mammalian Erythrocyte Micronucleus Test of MWCNT in ICR mice. Title

This test was performed to assess the ability of MWCNT to induce **Objectives**

chromosomal aberrations in CHO-k1 cells.

Sponsor Name : Bioconvergence Technology Laboratory,

Korea Conformity Laboratories

Address : 8, Gaetbeol-ro 145 beongil, Yeonsu-gu, Incheon, Korea

Tel. : +82-32-859-4046 Fax : +82-32-858-0020

Testing Name : Bioconvergence Technology Laboratory,

Korea Conformity Laboratories facility

> Address : 8, Gaetbeol-ro 145 beongil, Yeonsu-gu, Incheon, Korea

Tel. : +82-32-859-4046 Fax : +82-32-858-0020

Study Animal aquisition : 14 May 2013 **Schedule** Grouping : 20 May 2013 Administration : 21 May 2013 Slide preparation : 22 May 2013 Submission of final report : 23 2013

Archives

1) Period of storage: 5 years after study completion

2) Documents of storage: Study plan, Documents related to the test substance, Raw data, Final report, Documents related to the GLP

August

3) Room of specimen storage: Slides

4) Room of storage: CD, Related documents

Contents

1. Summary ·····	··· 1
2. Test substance and control substances ·····	2
3. Preparation of test substance	4
4. Test animal ·····	4
5. Breeding Environment	5
6. Test method ·····	6
7. Results ·····	· 10
1) Micronuclei frequency and cytotoxicity	· 10
2) Clinical sign ·····	· 10
3) Body weight ····	
8. Discussion and conclusion ·····	· 10
9. References ·····	· 11
10. TABLES (Group summary)	
Table 1. Test results (main study)	
Table 2. Body weights of animals (main study)	· 13
11. Appendices (Individual data)	
Appendix 1. Test results (main study)	· 14
Appendix 2. Body weights of animals (main study)	· 15
12. Annexes	
Annex 1. Test results (preliminary range-finding study)	· 16
Annex 2. Body weights of animals (preliminary range-finding study)	· 17
Annex 3. Test substance data sheet ····	
Annex 4. Certificate of drinking water analysis	· 19
Annex 5. Laboratory diet certification report	
Annex 6. Diagnostic summary report of test animals	
Annex 7. Certificate of animal room	· 28
Annex 8. GLP Certificate	
Annex 9. Quality assurance statement-Original ·····	
Annex 10. Study personnel-Original ·····	. 34

1. Summary

The test substance, MWCNT, was evaluated for its potential to induce micronucleus in bone marrow of the oral administrated mouse.

Animals were 8 weeks old at administration of test substance. Test animals were dosed once orally with 3 dosing levels (50 mg/kg, 25 mg/kg, 12.5 mg/kg) and then frequency of MNPCEs and cytotoxicity (PCE ratio in the erythrocytes) was determined at 24 hours after administration of test substance. As a results of counting the frequency of MNPCEs of 2,000 PCEs, there was not statistically significant increase at any dose groups compare with negative control group. Also, there was statistically significant increase in positive control group compare with negative control group(p<0.05). The PCEs ratio of 200 erythrocytes (PCE/(PCE+NCE)), as a index of cytotoxicity, was not significantly decrease in PCE/(PCE+NCE) ratio compared to negative control group.

It was concluded that, under the condition of this study, test substance MWCNT did not induce micronucleus in mice bone marrow.

2. Test substance and control substances

- 1) Test substance
- (1) Product name: MWCNT (KUMHO: K-Nanos-100P)
- (2) Lot No.: Not available
- (3) Received date: 25 January, 2013
- (4) Received quantity: 666.89 g (including container weight)
- (5) Appearance: Powder
- (6) Purity: More than 90 % (We assumed 100 % for test substance and carried out the test)
- (7) Solubility : The test substance was dispersed in DPPC solution at 1.0 % concentration.
- (8) Storage condition: Room temperature
- (9) Stability: Not available
- (10) Caution: Not available
- (11) Supplier: KUMHO PETROCHEMICAL
- 2) Negative control substance (solvent 1)
- (1) Name: 1,2-Dipalmitoyl-sn-glycero-3-phosphocholine (DPPC)
- (2) Lot No.: 078K5203
- (3) CAS No.: 63-89-8
- (4) Molecular weight: 734.04
- (5) Date received: 21 March 2012
- (6) Quantity Received: 1 g
- (7) Appearance: White powder
- (8) Grade : \geq 99 %
- (9) Storage condition: Freezing storage
- (10) Supplier: Sigma-Aldrich, Inc.
- 3) Negative control substance (solvent 2)
- (1) Name: Dulbecco's phosphate buffered saline (DPBS)
- (2) Lot No.: 031M8307
- (3) CAS No.: Not available
- (4) Date received: 08 May 2012
- (5) Quantity Received: 480 g
- (6) Appearance: White powder
- (7) Storage condition: Cold storage
- (8) Supplier: Sigma-Aldrich, Inc.
- 4) Negative control substance (solvent 3)

Name: D-(+)-Glucose
 Lot No.: 071M0145V
 CAS No.: 50-99-7

(4) Date received: 28 August 2012

(5) Quantity Received : 1 kg(6) Appearance : White powder

(7) Storage condition: Room temperature

(8) Supplier: Sigma-Aldrich, Inc.

5) Negative control substance (solvent 4)

(1) Name: Bovine serum albumin

(2) Lot No.: 750462

(3) CAS No.: Not available

(4) Date received: 06 April 2009

(5) Quantity Received: 100 g

(6) Appearance: Brown powder

(7) Storage condition: Cold storage

(8) Supplier: Gibco

6) Justification for the selection of solvent

Prior to execution of this study, the solubility of the test substance was not well-dispersed in general dispersion agents. For this reason, we selected the DPPC solution (5.5 mM D-(+)-glucose + 0.6 mg/ml Bovine serum albumin + 0.01 mg/kg DPPC in DPBS) according to Kim et al., 2011 study (Evaluation of biocompatible dispersants for carbon nanotube toxicity tests, Arch. Toxicol. 85: 1499-1508). As a result of solubility test, we observed that the test substance was well-dispersed in DPPC solution at 1.0 % concentration.

7) Positive control substance

(1) Name: Mitomycin C

(2) CAS No.: 50-07-7

(3) Manufacturer : Sigma

(4) Product No.: M-4287

(5) Lot No.: SLBB7481V

(6) Character: Water soluble

(7) Storage condition : 4 $^{\circ}$ C

(8) Justification for the selection

Positive control substances were selected according to OECD guidelines No. 474.

8) Storage and treatment (KCL/SAM/013)

During the study, test substances were stored in the storage room No. 108-2. All

storage and treatment conditions are in accordance with the relevant regulations. The stability test and the homogeneity test were not performed because the test solution was prepared fresh in treatment day.

3. Preparation of test substance

The test substance was dispersed in DPPC solution and serially diluted with DPPC solution. The administration level of prepared test substance is determined to 20 ml/kg based on body weight measured at the administration date. The stability test and the homogeneity test were not performed because the test solution was prepared fresh in treatment day.

4. Test animal

- 1) Species and strains: SPF (Specific Pathogenic Free), CrljOri:CD1(ICR) mice
- 2) Supplier: Orient Bio Co., Ltd. (Address; 143-1, Sangdaewon-dong, Jungwon-gu, Seongnam, Gyeonggi-do, Korea)
- 3) Reason for animal selection: ICR mice has been widely used in mammalian erythrocyte micronucleus test. It has been known for suitable experimental animal for general toxicity. Moreover, sufficient fundamental data have been accumulated in toxicological study using ICR mice. Such data will be useful for interpretation and evaluation of test results.
- 4) Date of animal acquisition

Preliminary dose-range finding study: 02 April 2013

Main study: 14 May 2013

5) Number of animal

Preliminary dose-range finding study: 33 male mice

Main study: 33 male mice

6) Age at the arrival: 7 weeks old

7) Quarantine and acclimation

Preliminary dose-range finding study: After 7 days of the acquisition

Main study: After 7 days of the acquisition

(The only healthy animals were selected for this study by observation of general symptoms during the acclimation period)

8) Age of test animal at administration: 8 weeks old

9) Number of test animal

Preliminary dose-range finding study: 30 male mice

Main study: 30 male mice

- 10) Grouping Method: Animals were weighted one day before the administration, and grouped by using the graded body weight randomly.
- 11) Identification of individual animals: Individual animals was identified by tail marking using oil marker and individual cages were identified by card labeling. The record sheet was posted at the entrance of animal room and the study number, test title, test period, name of study director was recorded.
- 12) Weighing method of test animals: Body weight of test animals was measured at the time of animal acquisition, grouping and administration of test substance. Body weight was measured one by one in accordance with individual animal number using electrical balance calibrated by calibration standard weight and recorded in data sheet.
- 13) Disposal of remaining animals: The remaining animals were sacrificed by mercy killing. They are treated by SOP of this testing facility.
- 14) Ethics for animal use: This study was approved by the Korea Conformity Laboratories (KCL) Institutional Animal Care and Use Committee (IA13-00255).

5. Breeding Environment

- 1) Room Number: SPF Room No. 4 in the SPF animal facility area
- 2) Range of acceptable temperature and humidity

Preliminary dose-range finding study : 21.8 \pm 0.7 $^{\circ}$ C of temperature, 44.9 \pm 5.1 % of relative humidity

Main study : 21.5 ± 1.1 °C of temperature, 53.8 ± 5.1 % of relative humidity

- 3) Ventilation frequency : $10 \sim 15$ times/hr.
- 4) Lighting cycle: 12 hours of lighting duration (lighting up at 8 a.m. ~ lighting out at 8 p.m.)
- 5) Luminous intensity

Preliminary dose-range finding study: 284 Lux

Main study: 284 Lux

- 6) Breeding cage (size): During the quarantine, acclimation, administration and observation period, all animals were housed in polycarbonate cage (170 W x 235 L x 125 H mm)
- 7) Animal number per cage

Preliminary dose-range finding study: All animals were housed 5 animals per cage during quarantine and acclimation period and 3 animals per cage during administration and observation period

Main study: All animals were housed 5 animals per cage during quarantine and

acclimation period and 3 animals per cage during administration and observation period 8) Feed, water supply and litter

Feed (Teklad Certified Irradiated Global 18 % Protein Rodent Diet, Harlan Co. Ltd., USA) was supplied from DooYeol Biotech co., Ltd. and were provided free ingestion using feeding supplier. Incheon, Korea municipal tap water purified by reverse osmosis filtering system was provided ad libitum using water polycarbonate bottles. Litter was supplied from DooYeol Biotech co., Ltd. and used after autoclaved.

6. Test method

- 1) Administration (KCL/AEP/001)
- (1) Route of administration

Test substance was administrated by Intraperitoneal (I.P.).

(2) Frequency and duration of administration

Single dose at the morning of administration day

(3) Calculation of dosing volume

Individual dosage is adjusted based on the fasted body weight measured right before the administration. Dosing volume is 20 ml/kg body weight.

- 2) Determination of dose level
- (1) Preliminary dose-range finding study

Preliminary dose-range finding study was conducted for determination of dose levels and time of harvest in main study. Test groups were consisted of one negative control group, three dosing groups (12.5 mg/kg, 25 mg/kg, 50 mg/kg) and one positive control group (Mitomycin C (MMC) 2.0 mg/kg). Each group was consisted of three mice, harvested at 24 and 48 hours after administration of test substance and conducted for acute toxicity test and cytotoxicity in bone marrow.

(2) Main test

According to determined dose in preliminary study, the highest dose in the main study was 2000 mg/kg. Test groups were consisted of one negative control group (corn oil only), three dosing groups (12.5 mg/kg, 25 mg/kg, 50 mg/kg) and one positive control group (Mitomycin C (MMC) 2.0 mg/kg). Harvest of bone marrow was conducted at 24 hours after administration of test substances. Each group was consisted of six mice.

Group	Sex	Number of	Identification of	dosing volume
Group	Sex	animal	Animal	(ml/kg)
G1 (V.C)	Male	6	1~6	20
G2	Male	6	7~12	20
G3	Male	6	13~18	20
G4	Male	6	$19 \sim 24$	20
G5 (P.C)	Male	6	$25 \sim 30$	10

G1(V.C): Vehicle Control, G2~G4: Test group, G5(P.C): Positive Control

3) Reason for determination of dosage

To determined dosage, preliminary study was performed. In the preliminary test, test substance was administrated in dose of 50 mg/kg, 25 mg/kg and 12.5 mg/kg by oral. As the result, There was not observed dead animal and special abnormality compared to negative control group. When observing according to dose level (12.5 mg/kg, 25 mg/kg, 50 mg/kg) and harvest time after administration (24 and 48 hours), no significant decrease in PCE/(PCE+NCE) ratio compared to negative control group was observed (Annex 1). On the basis of preliminary test result, the dose levels in the main study were determined at 50 mg/kg (maximum dose), 25 mg/kg (middle dose) and 12.5 mg/kg (minimum dose) and harvest time was conducted at 24 hours after administration of test substance.

4) Test method

(1) Extraction of bone marrow and preparation of slides

At the appropriate harvest time points, the thigh bone of the test animal was collected with care to avoid blood contamination in autopsy room after putting the cervical vertebra out. The bone marrow was collected in a centrifuge tube by flushing the thigh bone inside with FBS (Hyclone, Lot No. AXC36539, AXC47493) The extracted bone marrow was centrifuged at 1,000 rpm for 5 minutes and then resuspended with small aliquots of FBS after discarding the supernatant. The bone marrow was smeared on a slide glass and then dried at room temperature and fixed in methanol for 5 minutes. Three slides a bone marrow were prepared.

(2) Staining of specimen

For the scoring the PCEs to NCEs ratio, the slides were stained in 4 % Giemsa staining solution. For the scoring the MNPCEs from PCEs, the slides were fixed and then stained with acridine orange (40 μ g/m ℓ) and covered with a cover glasses.

(3) Slide analysis

① Observation method

The observation of slides was performed with blind method and PCE to NCE ratio was examined by optical microscope at 400× magnification. The observation of MNPCEs was performed by fluorescence microscope equipped with FITC filter.

② Criteria

The PCE/NCE ratio was determined by scoring the number of PCEs and NCEs observed while scoring 200 erythrocytes per animal. The micronuclei frequency (expressed as percent micronucleated cells) was determined by analyzing the number of MNPCEs from 2,000 PCEs per animal. In case of observation by fluorescence microscope, PCE was counted cells which appeared as red-fluorescence color without nucleus by acridine orange staining. NCE stained with acridine orange was appeared as only shadow entity without fluorescence. In case of observation by optical microscope, PCE stained with Giemsa staining solution was counted cells appeared as purple or blue color. NCE stained with Giemsa staining solution was counted cells appeared as pink color. According to criteia in micronucleus judgement, the largest size was defined as 1/2 size of erythrocyte diameter and the smallest size was defined as up to the limit of identification. Shape was Including circle, donut, semi-circle. Color was defined as same color as nucleus of cell.

(4) Observation of general symptoms

It was conducted once at the administration date and autopsy date for observation of dead and abnormal sign of test animals

- (5) Measurement of body weight
- It was measured at the administration date and autopsy date.
- 5) Statistical analysis
- (1) Micronuclei frequency and PCE/(PCE+NCE) ratio
- ① Comparison of the negative control and treated groups: ANOVA test
- ② In case of p < 0.05 at 1), test for dose-response: Linear logistic regression
- ③ The negative and positive control groups: ANOVA test
- (2) Body weight
- ① Body weight of test animal at necropsy: ANOVA test
- ② In case of significant difference at ①: Dunnett T3 or Duncan's multiple range test

6) Criteria

In case PCE/(PCE+NCE) ratio (Mean±SD) is more than 0.2 at all the animals, test is judged to be reasonable. Determined as a positive result if there is dose-related increase in the number of micronucleated cells or clear increase in the number of

micronucleated cells in the single dose group and statistical significant increase in frequency of MNPCE(MNPCE/2000PCEs, Mean±SD, %). The result of statistical evaluation was regarded as significant when the P value was less than 0.05. Statistical analysis was performed with SPSS 12.0 program.

7. Results

1) Micronuclei frequency and cytotoxicity (Table 1. and Appendix 1.)

The frequencies of MNPCEs in 2,000 PCEs per animal were 0.08 %, 0.09 %, 0.08 %, 0.04 % and 5.73 % in order of negative control group, 12.5, 25, 50 mg/kg administration group and positive control group. There was no statistically increase in frequency of PCEs having micronuclei in administration group compared to negative control group. There was clear statistically significant increase in positive control group compared with negative control group in frequency of MNPCE (p<0.05).

compared with negative control group in frequency of MNPCE (p < 0.05).

The PCE ratio of 200 erythrocytes (PCE+NCE), as a index of cytotoxicity, was 0.39, 0.36, 0.36, 0.36, 0.41 and 0.28 in order of negative control, 12.5, 25, 50 mg/kg group and positive control. There was no significant decrease in PCE/(PCE+NCE) ratio

compared to negative control group.

2) Clinical sign

As a result of observation, no special abnormality was observed compared to negative control group at necropsy date.

3) Body weight (Table 2. and Appendix 2.)

There was no statistically change in body weight at any group compared to negative control group.

8. Discussion and conclusion

The mammalian erythrocyte micronucleus test was performed to evaluate genotoxicity of the test substance MWCNT using male mouse (ICR mouse).

Animals (ICR mouse) were 7 weeks old at acquisition and acclimated for more than 5 days before treatment of test substance. Test substance was dispersed in distilled water. The dose levels administered to the animals were as follows:

Preliminary does-range finding study: 500, 1000, 2000 mg/kg

Main study: 12.5, 25, 50 mg/kg

In preliminary dose-range finding study, there was not observed death of test animal in all dose level (12.5, 25, 50 mg/kg). When observing according to the dose level and for the harvest times (24, 48 hours) after administration, PCE/(PCE+NCE) ratio was not observed significant decrease in all administrated groups (12.5, 25, 50 mg/kg) compared

with negative control group.

In main study, the highest dose was 2000 mg/kg and harvest time point was 24 hours respectively which was referred to result of preliminary dose-range finding study.

There was no statistically significant in frequency of 2,000 PCEs having micronuclei in all test substance groups compared to negative control group. On the other hand, positive control group appeared statistically significant increase compared with negative control group (p<0.05). Also, there was no significant decrease in the number of PCEs per PCE+NCE was observed at the all dose levels (12.5, 25, 50 mg/kg) administered to the animals as compared to negative control.

Statistically change in body weight was not observed in all test group compared to negative control group.

It was concluded that, under the condition of this study, test substance MWCNT did not induce micronucleus in mice bone marrow.

9. References

- 1) National Institute of Environment Research (NIER) [Notice No. 2012-23, (revised 22 August 2013)].
- 2) OECD Guidelines for the Testing of Chemicals No. 474 "Mammalian Erythrocyte Micronucleus test" (Adopted: 21st July 1997)
- 3) Schmid, W. (1975). The Micronucleus Test. Mutation Res. 31: 9-15.
- 4) 林 眞 (1991). 小核試驗. -實驗法から データの 評價まで 醫藥安全性研究會 monograph series No. 2. サイエンテイスト社

10. Tables (Group Summary)

Table 1. Test results (main study)

Sampling time(hrs)	Groups	Dose (mg/kg)	Animal No.		PCE/2000 PCEs Mean±SD, %)	PCI (E/(PCE+NCE) Mean±SD)
		0	1	0.1		0.34	
			2	0.1		0.40	
	Wahiala aantaal		3	0		0.41	
	Vehicle control		4	0.25		0.44	
			5	0		0.37	
			6	0	(0.08 ± 0.10)	0.36	(0.39 ± 0.04)
		12.5	7	0.15		0.34	
			8	0.05		0.33	
			9	0.1		0.35	
			10	0.1		0.37	
			11	0.05		0.41	
			12	0.1	(0.09 ± 0.04)	0.36	(0.36 ± 0.03)
		25	13	0.05		0.32	
			14	0.1		0.41	
24	MANICNIT		15	0.15		0.40	
24	MWCNT		16	0		0.33	
			17	0.05		0.34	
			18	0.1	(0.08 ± 0.05)	0.37	(0.36 ± 0.04)
		50	19	0		0.51	
			20	0.05		0.38	
			21	0		0.50	
			22	0		0.44	
			23	0.1		0.37	
			24	0.1	(0.04 ± 0.05)	0.27	(0.41 ± 0.09)
		2.0	25	5.15		0.33	
			26	5.25		0.29	
	MMC		27	4.9		0.28	
	MMC		28	7.2		0.26	
			29	7.85		0.24	
			30	4	(5.73 ± 1.48) *	0.27	(0.28 ± 0.03) *

Vehicle: DPPC solution

* Significantly different from the control at P<0.05 (One-way ANOVA)

Abbreviations

MNPCE: PCE with one or more micronuclei

PCE : Polychromatic erythrocyte NCE : Normochromatic erythrocyte

MMC : Mitomycin C

Table 2. Body weights of animals (main study)

Sampling time(hrs) Groups	Dose Animal (mg/kg) No.		Body weights (gram, Mean±S.D)		
			Administration	Sacrifice	
	Vehicle control	0	6	33.67±1.22	33.80±1.37
		12.5	6	33.51±1.08	32.64±1.05
24	MWCNT	25	6	33.55±0.90	32.43±0.79
		50	6	33.22±0.60	33.10±0.50
	MMC	2.0	6	33.08±0.76	33.26±0.66

Vehicle : DPPC solution MMC : Mitomycin C

11. Appendices (Individual data)

Appendix 1. Test results (main study)

Sampling time(hrs)	Groups	Dose (mg/kg)	Animal No.	Slide No.	MNPCEs	PCEs	No PCE/	of NCE	Ratio
		0	1	19	2	2000	67	133	0.34
			2	30	2	2000	80	120	0.40
	Vehicle control		3	5	0	2000	82	118	0.41
	venicle control		4	23	5	2000	88	112	0.44
			5	3	0	2000	74	126	0.37
			6	29	0	2000	72	128	0.36
		12.5	7	14	3	2000	68	132	0.34
			8	4	1	2000	66	134	0.33
			9	22	2	2000	70	130	0.35
			10	6	2	2000	74	126	0.37
			11	20	1	2000	82	118	0.41
			12	15	2	2000	71	129	0.36
		25	13	10	1	2000	64	136	0.32
			14	25	2	2000	81	119	0.41
24	MWCNT		15	21	3	2000	80	120	0.40
24	WIWCNI		16	1	0	2000	65	135	0.33
			17	18	1	2000	67	133	0.34
			18	13	2	2000	73	127	0.37
		50	19	26	0	2000	101	99	0.51
			20	2	1	2000	76	124	0.38
			21	12	0	2000	100	100	0.50
			22	17	0	2000	87	123	0.44
			23	7	2	2000	74	126	0.37
			24	28	2	2000	54	146	0.27
		2.0	25	11	103	2000	65	135	0.33
			26	16	105	2000	58	142	0.29
	MMC		27	24	98	2000	56	144	0.28
	MMC		28	8	144	2000	52	148	0.26
			29	27	157	2000	47	153	0.24
			30	9	80	2000	54	146	0.27

Vehicle: DPPC solution

Abbreviations

MNPCE: PCE with one or more micronuclei

PCE : Polychromatic erythrocyte NCE : Normochromatic erythrocyte

MMC: Mitomycin C

Appendix 2. Body weights of animals (main study)

time(hrs) Chours (mg/kg) No. Administration Sacrifice 1 31.65 31.38 2 34.15 34.61 Vehicle control 3 33.18 33.37 5 34.39 34.35 6 35.19 33.67±1.22 35.35 33.80±1.37 7 32.81 32.14 8 32.33 31.84 33.25 10 34.08 33.25 11 34.20 32.96 12 35.03 33.51±1.08 34.25 32.64±1.05 13 31.95 31.43 14 33.51 32.25 32.64±1.05 16 33.91 33.71 32.25 16 33.91 33.71 32.43±0.79 24 MWCNT 25 16 33.91 33.71 17 33.68 32.18 32.18 18 34.72 33.55±0.90 32.91 32.43±0.79 20 32.36	Sampling	Groups	Dose	Animal		Body weigh	nt (gram)	
Vehicle control 0	time(hrs)	Groups	(mg/kg)	No.	Admi	nistration	Sac	crifice
Vehicle control 0 3 33.18 33.37 5 34.39 34.35 6 35.19 33.67±1.22 35.35 33.80±1.37 7 32.81 32.14 8 32.33 31.84 12.5 9 32.62 31.39 10 34.08 33.25 11 34.20 32.96 12 35.03 33.51±1.08 34.25 32.64±1.05 13 31.95 31.43 14 33.51 32.25 15 33.50 32.08 16 33.91 33.71 17 33.68 32.18 18 34.72 33.55±0.90 32.91 32.43±0.79 19 32.64 32.45 20 32.36 32.50 33.39 33.35 33.17 30 32.25 33.36 33.42 23 33.89 33.56 32.40 24 33.72 33.2±0.60 33.48 33.10±0.50 25 32.47 32				1	31.65		31.38	
Control 4 33.48 33.71 5 34.39 34.35 6 35.19 33.67±1.22 35.35 33.80±1.37 7 32.81 32.14 8 32.33 31.84 12.5 9 32.62 31.39 10 34.08 33.25 11 34.20 32.96 12 35.03 33.51±1.08 34.25 32.64±1.05 13 31.95 31.43 14 33.51 32.25 15 33.50 32.08 16 33.91 33.71 17 33.68 32.18 18 34.72 33.55±0.90 32.91 32.43±0.79 19 32.64 32.45 20 32.36 32.50 21 33.35 33.17 22 33.36 33.42 23 33.89 33.56 24 33.72 33.22±0.60 33.48 33.10±0.50 MMC 2 28 33.35 33.60 29 32.75 33.10				2	34.15		34.61	
Control 4 33.48 33.71 5 34.39 34.35 6 35.19 33.67±1.22 35.35 33.80±1.37 7 32.81 32.14 8 32.33 31.84 12.5 9 32.62 31.39 10 34.08 33.25 11 34.20 32.96 12 35.03 33.51±1.08 34.25 32.64±1.05 13 31.95 31.43 14 33.51 32.25 15 33.50 32.08 16 33.91 33.71 17 33.68 32.18 18 34.72 33.55±0.90 32.91 32.43±0.79 19 32.64 32.45 20 32.36 32.50 21 33.35 33.17 22 33.36 32.50 50 21 33.35 33.17 22 33.36 33.42 23 33.89 33.56 24 33.72 33.22±0.60 33.48 33.10±0.50 MMC 2 28 33.35 33.60 29 32.75 33.19		Vehicle	0	3	33.18		33.37	
6 35.19 33.67±1.22 35.35 33.80±1.37 7 32.81 32.14 8 32.33 31.84 12.5 9 32.62 31.39 10 34.08 33.25 11 34.20 32.96 12 35.03 33.51±1.08 34.25 32.64±1.05 13 31.95 31.43 14 33.51 32.25 15 33.50 32.08 16 33.91 33.71 17 33.68 32.18 18 34.72 33.55±0.90 32.91 32.43±0.79 19 32.64 32.45 20 32.36 32.50 21 33.35 33.17 22 33.36 32.50 21 33.35 33.17 22 33.36 33.42 23 33.89 33.56 24 33.72 33.22±0.60 33.48 33.10±0.50 MMC 2 28 33.35 33.60 29 32.75 33.19		control	U	4	33.48		33.71	
7 32.81 32.14 8 32.33 31.84 12.5 9 32.62 31.39 10 34.08 33.25 11 34.20 32.96 12 35.03 33.51±1.08 34.25 32.64±1.05 13 31.95 31.43 14 33.51 32.25 15 33.50 32.08 16 33.91 33.71 17 33.68 32.18 18 34.72 33.55±0.90 32.91 32.43±0.79 19 32.64 32.45 20 32.36 32.50 21 33.35 33.17 22 33.36 32.50 21 33.35 33.17 22 33.36 33.42 23 33.89 33.56 24 33.72 33.22±0.60 33.48 33.10±0.50 MMC 2 26 32.42 32.59 27 33.06 32.97 28 33.35 33.60 29 32.75 33.19				5	34.39		34.35	
8 32.33 31.84 12.5 9 32.62 31.39 10 34.08 33.25 11 34.20 32.96 12 35.03 33.51±1.08 34.25 32.64±1.05 13 31.95 31.43 14 33.51 32.25 15 33.50 32.08 16 33.91 33.71 17 33.68 32.18 18 34.72 33.55±0.90 32.91 32.43±0.79 19 32.64 32.45 20 32.36 32.50 21 33.35 33.17 22 33.36 32.50 21 33.35 33.17 22 33.36 33.42 23 33.89 33.56 24 33.72 33.22±0.60 33.48 33.10±0.50 25 32.47 32.79 26 32.42 32.59 MMC 2 28 33.35 33.60 29 32.75 33.19				6	35.19	33.67±1.22	35.35	33.80±1.37
12.5 9 32.62 31.39 10 34.08 33.25 11 34.20 32.96 12 35.03 33.51±1.08 34.25 32.64±1.05 13 31.95 31.43 14 33.51 32.25 15 33.50 32.08 16 33.91 33.71 17 33.68 32.18 18 34.72 33.55±0.90 32.91 32.43±0.79 19 32.64 32.45 20 32.36 32.50 31.43 50 21 33.35 33.17 22 33.36 32.50 24 33.72 33.22±0.60 33.48 33.10±0.50 25 32.47 32.79 26 32.42 32.59 27 33.06 32.97 MMC 2 28 33.35 33.60 29 32.75 33.19				7	32.81		32.14	
10 34.08 33.25 11 34.20 32.96 12 35.03 33.51±1.08 34.25 32.64±1.05 13 31.95 31.43 14 33.51 32.25 15 33.50 32.08 16 33.91 33.71 17 33.68 32.18 18 34.72 33.55±0.90 32.91 32.43±0.79 19 32.64 32.45 20 32.36 32.50 21 33.35 33.17 22 33.36 33.42 23 33.89 33.56 24 33.72 33.22±0.60 33.48 33.10±0.50 25 32.47 32.79 26 32.42 32.59 MMC 2 28 33.35 33.60 29 32.75 33.19				8	32.33		31.84	
11 34.20 32.96 12 35.03 33.51±1.08 34.25 32.64±1.05 13 31.95 31.43 14 33.51 32.25 24 MWCNT 25 15 33.50 32.08 16 33.91 33.71 17 33.68 32.18 18 34.72 33.55±0.90 32.91 32.43±0.79 19 32.64 32.45 20 32.36 32.50 21 33.35 33.17 22 33.36 33.42 23 33.89 33.56 24 33.72 33.22±0.60 33.48 33.10±0.50 25 32.47 32.79 26 32.42 32.59 27 33.06 32.97 28 33.35 33.60 29 32.75 33.19			12.5	9	32.62		31.39	
12 35.03 33.51±1.08 34.25 32.64±1.05 13 31.95 31.43 14 33.51 32.25 15 33.50 32.08 16 33.91 33.71 17 33.68 32.18 18 34.72 33.55±0.90 32.91 32.43±0.79 19 32.64 32.45 20 32.36 32.50 21 33.35 33.17 22 33.36 33.42 23 33.89 33.56 24 33.72 33.22±0.60 33.48 33.10±0.50 25 32.47 32.79 26 32.42 32.59 MMC 2 28 33.35 33.60 29 32.75 33.19				10	34.08		33.25	
13 31.95 31.43 14 33.51 32.25 15 33.50 32.08 16 33.91 33.71 17 33.68 32.18 18 34.72 33.55±0.90 32.91 32.43±0.79 19 32.64 32.45 20 32.36 32.50 21 33.35 33.17 22 33.36 33.42 23 33.89 33.56 24 33.72 33.22±0.60 33.48 33.10±0.50 25 32.47 32.79 26 32.42 32.59 27 33.06 32.97 28 33.35 33.60 29 32.75 33.19			11	34.20		32.96		
14 33.51 32.25 15 33.50 32.08 16 33.91 33.71 17 33.68 32.18 18 34.72 33.55±0.90 32.91 32.43±0.79 19 32.64 32.45 20 32.36 32.50 21 33.35 33.17 22 33.36 33.42 23 33.89 33.56 24 33.72 33.22±0.60 33.48 33.10±0.50 25 32.47 32.79 26 32.42 32.59 MMC 2 28 33.35 33.60 29 32.75 33.19				12	35.03	33.51 ± 1.08	34.25	32.64±1.05
MWCNT 25				13	31.95		31.43	
MWCNT 25 16			NT 25	14	33.51		32.25	
16 33.91 33.71 17 33.68 32.18 18 34.72 33.55±0.90 32.91 32.43±0.79 19 32.64 32.45 20 32.36 32.50 21 33.35 33.17 22 33.36 33.42 23 33.89 33.56 24 33.72 33.22±0.60 33.48 33.10±0.50 25 32.47 32.79 26 32.42 32.59 MMC 2 28 33.35 33.60 29 32.75 33.19	24	MWCNT		15	33.50		32.08	
$ \begin{array}{c ccccccccccccccccccccccccccccccccccc$		IVIWCNI		16	33.91		33.71	
19 32.64 32.45 20 32.36 32.50 21 33.35 33.17 22 33.36 33.42 23 33.89 33.56 24 33.72 33.22±0.60 33.48 33.10±0.50 25 32.47 32.79 26 32.42 32.59 27 33.06 32.97 28 33.35 33.60 29 32.75 33.19				17	33.68		32.18	
20 32.36 32.50 21 33.35 33.17 22 33.36 33.42 23 33.89 33.56 24 33.72 33.22±0.60 33.48 33.10±0.50 25 32.47 32.79 26 32.42 32.59 27 33.06 32.97 28 33.35 33.60 29 32.75 33.19				18	34.72	33.55±0.90	32.91	32.43±0.79
50 21 33.35 33.17 22 33.36 33.42 23 33.89 33.56 24 33.72 33.22±0.60 33.48 33.10±0.50 25 32.47 32.79 26 32.42 32.59 27 33.06 32.97 28 33.35 33.60 29 32.75 33.19				19	32.64		32.45	
MMC 2 33.36 33.42 33.56 24 33.72 33.22±0.60 33.48 33.10±0.50 25 32.47 32.79 26 32.42 32.59 27 33.06 32.97 28 33.35 33.60 29 32.75 33.19				20	32.36		32.50	
MMC 2 33.36 33.42 23 33.89 33.56 24 33.72 33.22±0.60 33.48 33.10±0.50 25 32.47 32.79 26 32.42 32.59 27 33.06 32.97 28 33.35 33.60 29 32.75 33.19			50	21	33.35		33.17	
MMC 2 24 33.72 33.22±0.60 33.48 33.10±0.50 25 32.47 32.79 26 32.42 32.59 27 33.06 32.97 28 33.35 33.60 29 32.75 33.19			30	22	33.36		33.42	
MMC 2 32.47 32.79 26 32.42 32.59 27 33.06 32.97 28 33.35 33.60 29 32.75 33.19				23	33.89		33.56	
MMC 2				24	33.72	33.22 ± 0.60	33.48	33.10±0.50
MMC 2 27 33.06 32.97 28 33.35 33.60 29 32.75 33.19				25	32.47		32.79	
MMC 2 28 33.35 33.60 29 32.75 33.19	N046			26	32.42		32.59	
28 33.35 33.60 29 32.75 33.19		MAC	2	27	33.06		32.97	
		MMC	2	28	33.35		33.60	
30 34.45 33.08±0.76 34.41 33.26±0.66				29	32.75		33.19	
				30	34.45	33.08±0.76	34.41	33.26±0.66

Vehicle : DPPC solution MMC : Mitomycin C

12. Annexes

Annex 1. Test results (preliminary range-finding study)

Sampling	Canara	Dose	Animal	MNPCE	/2000 PCEs	PCE/(P	CE+NCE)	
time(hrs)	Groups	(mg/kg)	(mg/kg) No. (Mean±SD,%)		n±SD,%)	(Mean±SD)		
	Vehicle		1	0.1		0.45		
	control	0	2	0.15		0.41		
	control		3	0.25	0.17 ± 0.08	0.39	0.42 ± 0.03	
			4	0.1		0.46		
		12.5	5	0.1		0.45		
			6	0.15	0.12 ± 0.03	0.42	0.44 ± 0.02	
			7	0.05		0.31		
24	MWCNT	25	8	0.05		0.44		
			9	0.15	0.08 ± 0.06	0.49	0.41 ± 0.10	
			10	0.1		0.45		
		50	11	0		0.45		
			12	0.1	0.07 ± 0.06	0.32	0.40 ± 0.08	
			13	5.25		0.42		
	MMC	2	14	5.6		0.45		
			15	5.55	5.47±0.19*	0.50	0.45 ± 0.04	
	Vehicle		16	0.1		0.41		
	control	0	17	0.25		0.47		
	Control		18	0.25	0.20±0.09	0.42	0.43 ± 0.03	
			19	0.2		0.48		
		12.5	20	0.1		0.51		
			21	0.05	0.12 ± 0.08	0.47	0.49 ± 0.02	
			22	0		0.54		
48	MWCNT	25	23	0.25		0.48		
			24	0.1	0.12 ± 0.13	0.38	0.47 ± 0.08	
			25	0.05		0.46		
		50	26	0.1		0.38		
			27	0.05	0.07±0.03	0.47	0.43±0.05	
			28	3.65		0.46		
	MMC	2	29	3.6		0.25		
			30	5.6	4.28±1.14*	0.26	0.32±0.12	

Vehicle: DPPC solution

* Significantly different from the control at P<0.05 (One-way ANOVA)

Abbreviations

MNPCE: PCE with one or more micronuclei

PCE : Polychromatic erythrocyte NCE : Normochromatic erythrocyte

MMC : Mitomycin C

Annex 2. Body weights of animals (preliminary range-finding study)

Sampling	Groups	Dose	Animal		Body weigh	nt (gram)		
time(hrs)	Groups	(mg/kg)	No.	Admi	nistration	Sac	Sacrifice	
	Vehicle		1	31.82		32.28		
	control	0	2	33.26		33.01		
	Control		3	34.06	33.05 ± 1.14	35.00	33.43±1.41	
			4	31.05		31.08		
		12.5	5	33.41		33.09		
			6	33.26	32.57 ± 1.32	32.87	32.35±1.10	
			7	32.92		33.17		
24	MWCNT	25	8	32.95		31.65		
			9	33.78	33.22 ± 0.49	34.11	32.98±1.24	
			10	32.36		30.60		
		50	11	32.54		30.87		
			12	34.72	33.21±1.31	33.08	31.52±1.36	
			13	32.24		32.27		
	MMC	2	14	32.71		33.05		
			15	33.74	32.90 ± 0.77	34.14	33.15±0.94	
	Vehicle		16	32.53		32.80		
	control	0	17	32.68		33.29		
	Control		18	34.65	33.29±1.18	33.81	33.30±0.51	
			19	31.73		31.79		
		12.5	20	32.43		32.78		
			21	34.11	32.76 ± 1.22	33.95	32.84±1.08	
			22	32.60		32.44		
48	MWCNT	25	23	32.00		32.38		
			24	33.93	32.84 ± 0.99	33.65	32.82±0.72	
			25	32.15		31.17		
		50	26	32.60		32.09		
			27	33.70	32.82 ± 0.80	33.33	32.20±1.08	
			28	31.96		32.69		
	MMC	2	29	32.22		32.44		
			30	35.09	33.09±1.74	35.74	33.62±1.84	

Vehicle : DPPC solution MMC : Mitomycin C

Annex 3. Test substance data sheet



Annex 4. Certificate of drinking water analysis

TEST REPORT

(1) 1.30
.30
Coe

TEST REPORT

No : PC12-00576

Test Items	Units	Limitations	LOQ	Test Results
Total colony counts	CFU/mL	Less than 100	0	34
Total coliforms	-/(100mL)	Not detected		Not detected
E-Coli	-/(100mL)	Not detected	-	Not detected
Lead	mg/L	Less than 0.01	0,005	Not detected
Arsenic	mg/L	Less than 0.01	0.005	Not detected
Selenium	mg/L	Less than 0.01	0.005	Not detected
Cadmium	mg/L	Less than 0.005	0.002	Not detected
Bolon	mg/L	Less than 1.0	0.01	Not detected
Copper	mg/L	Less than 1.0	0.008	Not detected
Zinc	mg/L	Less than 3.0	0.002	0,003
Iron	mg/L	Less than 0.3	0.05	Not detected
Manganese	mg/L	Less than 0.3	0.005	Not detected
Aluminium	mg/L	Less than 0.2	0.02	Not detected
Mercury	mg/L	Less than 0.001	0.001	Not detected
Fluoride	mg/L	Less than 1.5	0.15	Not detected
Nitrate nitrogen	mg/L	Less than 10	0.1	0.2
Chloride	mg/L	Less than 250	0.4	Not detected
Sulfate	mg/L	Less than 200	2	Not detected
Diazinon	mg/L	Less than 0.02	0.0005	Not detected
Parathion	mg/L	Less than 0.06	0.0005	Not detected
Femitrotion	mg/L	Less than 0.04	0.0005	Not detected
Dichloromethane	mg/L	Less than 0.02	0.002	Not detected
1.1.1-Trichloroethane	mg/L	Less than 0.1	0.001	Not detected
Benzene	mg/L	Less than 0.01	0.001	Not detected
Toluene	mg/L	Less than 0.7	0.001	Not detected
Ethylbenzene	mg/L	Less than 0.3	0.001	Not detected
Xylene	mg/L	Less than 0.5	0.001	Not detected
1,1-Dichloroethylene	mg/L	Less than 0.03	0.001	Not detected
Tetrachlorocarbon	mg/L	Less than 0.002	0.001	Not detected
Tetrachloroethylene	mg/L	Less than 0.01	0.001	Not detected
Trichloroethylene	mg/L	Less than 0.03	0.001	Not detected
,2-Dibromo-3-Chloropropane	mg/L	Less than 0.003	0.001	Not detected

- page 2 of 3

QP-20-01-08(1)

No : PC12-00576

TEST REPORT

Test Items	Units	Limitations	LOQ	Test Results
Carbaryl	mg/L	Less than 0.07	0.005	Not detected
Chromium	mg/L	Less than 0.05	0.003	Not detected
Ammonia nitrogen	mg/L	Less than 0.5	0.01	Not detected
PhenoI	mg/L	Less than 0.005	0.005	Not detected
Detergent	mg/L	Less than 0.5	0.1	Not detected
Cyani de	mg/L	Less than 0.01	0.01	Not detected
pH	1+	5.8 ~ 8.5	H)	6.4
Turbidity	NTU	Less than 1	0.02	0.16
Color	degree	Less than 5	1	Not detected
Taste	1.7-	Free	0	Pass
Odor	0.0	Free	-	Pass
Hardness	mg/L	Less than 300	1	Not detected
Consumption of KMnO4	mg/L	Less than 10	0.3	0.9
Total solids	mg/L	Less than 500	2	4
Test method	Notifical	tion No.2012-143 o	f the Mini	stry of Environ

--- End of Report ---

- page 3 of 3 -

QP-20-01-08(1)

Reissuance (R1)

Date : 2013.4.11

TEST REPORT

1. No: PC13-00284

2. Client

○ Name : Korea Conformity Laboratories(Incheon)

Address: #7-44, Songdo-dong, Yeonsu-gu, Incheon, Korea

Date of Receipt : Mar. 14, 2013
○ Date of Issued : Apr. 17, 2013
3. Use of Report : Submission

4. Test Sample : Drinking Water (Arrimal room)

5. Method :

(1) Notification No.2012-143 of the Ministry

of Environment

Aftirmation Tested By Name : Hyoung jun Seok Seok. Technical Manager Name : Sang Cheul Lee S. C. Ge

The case of acrols only to the stitutents of processing identified and to the sample(s) tested unless otherwise specified.

The hard results are not indicative of representative of the qualifies of the capitifies of the lot from which the sample was taken or of apparently identified to which products.

Korea Conformity Laboratories

President Song Jae Bin

Jae Bin Song

Address : 704-932 277-5. Jukleon-Dong, Dalseo-Gui, Daegui, 704-932, Korea 82-53-557-6681

Result Inquiry: Environmental Testing Center 82-2-2102-2598

- Page 1 of 3

QP-20-01-07(2)

TEST REPORT

No. 1 PC13-00284

6 Test Results

1) Drinking Water (Arimal toom)

Test Item(s)	Unit	Limitation(s)	L00	Test method used	Test Result(a
Total colony counts	CFU/nL	Less than 100	0	(1)	0:
Total colliforms	-/(100mL)	Not detected	-	(1)	Not detected
E-Coll	-/(100mL)	Not detected	= -	(1)	Not detected
Lead	mg/L	Less than 0.01	0.005	(1)	Not detected
Arsenic	mg/L	Less than 0.01	0.005	(1)	Not detected
Selenium	mg/L	Less than 0.01	0.005	(1)	Not detected
Cadmilum	mg/L	Less than 0,005	0.002	(1)	Not detected
Boton	mg/L	Less than 1.0	0.01	(1)	Not detected
Copper	mg/L	Less than 1.0	800.0	(1)	Not detected
Zinc	mg/L	Less than 3.0	0.002	(1)	0.003
Lron	mg/L	Less than 0.3	0.05	(1)	Not detected
Manganese	mg/L	Less than 0.3	0.005	(1)	Not detected
Alum(nlum	mg/L	Less than 0.2	0.02	(1)	Not detected
Mercury	mg/L	Less than 0.001	0:001	(1)	Not detected
Fluoride	mg/L	Less than 1.5	0.15	(1)	Not detected
Nifrate nitrogen	rig/L	Less than 10	0.1	(1)	0.2
Chibride	mg/L	Less than 250	0.4	(1)	0.6
Sulfate	mg/L	Less than 200	2	(1)	Not detected
Diazinon	mg/L	Less than 0.02	0.0005	(1)	Not detected
Parathion	mg/L	Less than 0.06	0.0005	111	Not detected
Fenitrotion	mg/L	Less than 0.04	0.000.5	(1)	Not detected
Dichloromethane	mg/L	Less than 0.02	0.002	(1)	Not detected
1.1.1-Trichtoroethane	mg/L	Less than 0.1	0.001	(1)	Not detected
Benzene	mg/L	Less than 0.01	0.001	(1)	Not detected
Toluene	mg/L	Less than 0.7	0.001	(1)	Not detected
Ethylbenzene	mg/t	Less than 0.3	0.001	(1)	Not detected
XyTene	mg/L	Less than 0.5	0,001	(1)	Not detected
1.1-Dichloroethylene	mg/L	Less than 0.03	0.001	(1)	Not detected
Tetrachiorocarbon	mg/L	Less than 0.002	0.001	(1)	Not detected
Tetrach/orgethy/ene	mg/L	Less than 0.01	0.001	(1)	Not detected

Fage 2 at 3

QF-20-01-08(2)

TEST REPORT

No : PC13-00284

6. Test Results

1) Drinking Water (Animal room)

Test Item(s)	Unit	Limitation(s)	L00	Test method used	Test Result(s)
Trichloroethylene	mg/L	Less than 0.03	0.001	(1)	Not detected
1.2-Dibrono-3-Chloropropane	_mg/L	Less than 0.003	0.001	(1)	Not detected
Carbaryl	mg/L	Less than 0.07	0:005	(1)	Not detected
Chromium	mg/L	Less than 0.05	0.03	(1)	Not detected
Ammonium Nitrogen	mg/L	Less than 0:5	0.01	(1)	Not detected
Phenol	mg/L	Less than 0.005	0.005	(1)	Not detected
Alkyl Benzene Sulfate	mg/L	Less than 0.5	0.1	(1)	Not detected
Cyanide	mg/L	Less than 0.01	0.01	(1)	Not detected
pH	-	5.8 - 8.5	-	(1)	6.2
Turbidity	NTU	Less than 1	0.02	(1)	0.11
Color	degree	Less than 5	1	(1)	Not detected
Taste	9 1	Free	3	(1)	Pass
.Odor		Free	- 1	(1)	Pass
Hardness	mg/L	Less than 300	1	(1)	Not detected
Consumption of KMnO ₄	mg/L	Less than 10	0.3	(4)	0.6
Total solids	mg/L	Less than 500	2	(1)	Not detected

- End of Report -

Page 3 of 3

OP-20H0.1H08(2)

Annex 5. Laboratory diet certification report

Laboratory Diet Certification Report

Teklad Certified Irradiated Global 18% Protein Rodent Diet Lot Number 2918C-120212MA 2918C Date of Manufacture 12/02/12 Report Date 12/18/12 The following data is a consolidation of results obtained **Proximate Analysis** from one or more independent testing laboratories. The Recutt (%) actual laboratory results are available upon request. Protein 18.40 Fat 6.14 I have reviewed this Fiber 3.32 Kut Schafer document Moisture 12.00 Ash 2012.12.21 07:18:28 5.51 Calcium 0.98 -06'00' Phosphorus 83.0 Feed Contaminant Screen Established Maximum Analysis Heavy Metals 0.17 Arsenic 1.00 Cadmium < 0.10 0.50 ppm < 0.20 Lead 1.50 ppm Mercury < 0.05 ppm 0.20 Selenium 0.24 0.50 ppm Mycotoxin Aflatoxin B1, B2, G1, G2 < 5.00 5.00 ppb Chlorinated Hydrocarbons Aldrin < 0.01 0.03 < 0.01 Lindane 0.05 ppm Chlordane < 0.01 0.05 ppm DDT & related substances < 0.03 0.15 ppm Dieldrin < 0.02 0.03 ppm Endrin < 0.02 0.03 ppm Heptachion < 0.01 0.03 ppm Heptachlor Epoxide < 0.01 ppm 0.03 Toxaphene < 0.10 0.15 ppm PCB's < 0.10ppm 0.15 a-BHC < 0.01 0.05 ppm ь-внс < 0.01 ppm 0.05 d-BHC < 0.01 0.05 ppm Hexach orobenzene < 0.01 0.03 ppm Mirex < 0.01 ppm 0.02 Methoxychlor < 0.05 0.50 ppm Organophosphates < 0.15 0.50 Thimet ppm Diazinon < 0.14 0.50 Disulfaton < 0.15 0.50 ppm Methyl Parathion < 0.14 0.50 ppm Malathion < 0.14 0.50 Parathion < 0.12 0.50 ppm Thiodan < 0.02 ppm 0.50 Ethion < 0.14 0.50 ppm Trithion < 0.15 ppm 0.50 P.O. Box 44220 - Madison, WI 53744-4220 - 800 -483 -5523 - www.harlan.com harlan* Harles, Harles Laboratories, Hoping you to receive better, and the Harles logs are trademarks and trade sames of Harles Laboratories, by: 6 2006

Laboratory Diet Certification Report

Teklad Certified Irradiated Global 18% Protein Rodent Diet Lot Number 2918C-030413MA 2918C Date of Manufacture 03/04/13 Report Date 03/19/13 The following data is a consolidation of results obtained Proximate Analysia from one or more independent testing laboratories. The Recult (%) actual laboratory results are available upon request. Protein Fat 6.17 I have reviewed this Fiber 3.82 Kut Schofe document Moisture 10.50 Ash 5.66 2013.03.20 09:52:32 Quality Assurance Coordinate: Tekind Dis Research Models and Services Calcium 1.01 -05'00' Phosphorus 0.77 Feed Contaminant Screen Established Maximum Reput Concentration Heavy Metals Arsenic 0.12 1.00 < 0.10 Cadmium 0.50 ppm Lead < 0.20 1.50 ppm Mercury < 0.05 0.20 ppm 0.34 Selenium ppm 0.50 Mycotoxin Aflatoxin B1, B2, G1, G2 < 5.00 ppp 5.00 Chlorinated Hydrocarbons < 0.01 0.03 ppm Lindane < 0.01 0.05 Chlordane < 0.01 0.05 ppm DDT & related substances < 0.03 ppm 0.15 Dieldrin < 0.02 0.03 ppm Endrin < 0.02 0.03 ppm Heptachlor < 0.01 0.03 ppm Heptachlor Epoxide < 0.01 0.03 ppm Toxaphene < 0.10 0.15 ppm PCB's < 0.10 0.15 ppm a-BHC < 0.01 ppm 0.05 ь-внс < 0.01 0.05 ppm d-BHC < 0.01 ppm 0.05 Hexachlorobenzene < 0.01 0.03 ppm < 0.01 ppm 0.02 Methoxychlor < 0.05 ppm Organophosphates < 0.15 0.50 Diazinon < 0.14 0.50 ppm Disulfaton < 0.15 0.50 ppm Methyl Parathion < 0.14 0.50 ppm Malathion < 0.14 0.50 ppm

P.O. Box 44220. - Madison, WI 53744-4220 800 483 5523 - www.harlan.com

ppm

ppm

ppm

ppm

< 0.12

< 0.02

< 0.14

< 0.15

Parathion

Thiodan

Trithion

Ethion



0.50

0.50

0.50

0.50

CR Rodent Production

Sponsor: Orient Bio Inc.

Past 18 Months

Positive / Tesset

67.48

07.48

0 / 48


0 / 48

07.48

0 / 48

0 / 48.

Annex 6. Diagnostic summary report of test animals



REO se	MFIA.	2013-05	0/	8	0 /	48
(DIM et	MFIA.	2013-05	0.0	8	0.7	48
LCMV at	MEIA	2013-03	(17)	N .	0/	48
ECTRO as	MFIA	2013-05	0.0	8	0/	48
MHT we	META	2013-05	.07	N	07	48
MCMF as	MFIA.	2013-05	0.7	8	-07	48
X 40	MFIA	2011-05	-0/	N	0.7	48
POLT #	MEIA	2013-05	07	8	0/	48
HANT as	MEIA	201.1-05	0.7		0/	4N
MTT.V as	MEIA	2013-05	-01/	N	0.7	48
LDV at	MFIA	2012-41	0/	8	0/	16
	Microbiol	logy				
II bremekisepusa be	Culture	2013-05	97	X	0/	AN
CAR Bacillus as	MFIA/PCR	2013-03	01	8	0/	48
C rodenium ae	Culture	2013-05	-97	8.	07	36

Trezzer's Distance og	Exam	2013-05	0/	8	(0)	48
Strep pronounement be	Culture	2013-03	0.7	8	.0/	48.
S. mondeforms at	PCR	2013-05	0.7	8	0.0	12
Shiph attirous re-	Culture	2013-05	07	8		48
Salmonetin spp 🐠	Culture	2013-05	07	8	0.7	48
A surragement on	Culture	2013-05	07	8.	1.007	48
P ¹ , presumotropica oc	Culture	2013-05	0.7	8	07	48
P. multocida se	Culture	2013-05	07	8	0.7	48
NY pulminus ad	MFIA	2013-85	0.7	×	0/	48
K presimoniae se	Culture	2013-01	07	8	0.7	48
K incellucia aw	Culture	2013-05	0/		0.7	48
Helicobactor sp. be	PCR.	2013-05	- 0/	B	0.7	48
R. hepanicus - se	PCR.	2013-05	-07			48
D. Allia te	PCR.	2013-05	07	В		48
C hinscheri ae	Culture	2013-05	97	8	0/	48
C rodentium ne	Culture	2013-05	-07		0.6	36
CAR Bacillus as	MFIA/PCR	2013-03	07	8	.0./	48
II bromckisegment be	Culture	2013-05	97	X	0.7	AN

Vizzar's Districte ag	Exam	2013-05	0 / 8	0 / 48
	Pathology			
Grown Carm. +h	Exam, Histopathology	2013-05	0/-8	0 / 48
	Parasitology			
Ectoparasites =	Exam	2013-05	07.8	07.48
Helminika w	Exam	2013-05	07.8	07.48
Giardia sp. ne	Esam	2013-05	0.8	0 / 28
Spacemaclicus up he	Exam	2913-05	0.8	0.0.28
Officer Provinsian ee	Exam	2013-05	.07.8	0./ 28
E cuntrall to	MEIA	2013-05	0.8	0.5.48

I THEONY POLICY FOR POSITIVE RESULT: a - immediate termination; b - planted throse recycle of the colony, c - no action. TEST [NG SET REDICTE of summer development of the control parties of

Charles River RADS-ILIMS Report

Page I of I

Annex 7. Certificate of animal room

Certific	cation of Environment for animal breeding room
Study No.	GT13-00019
Title	Mammalian Erythrocyte Micronucleus Test of MWCNT in ICR mice (Preliminary test)
SPF Room No.	SPF #4 Animal Room
Period of animal Breeding	2013-04-02 ~ 2014-04-11

Breeding environment condition

Section	Range of SOP	Survey value	Remark
Temperature	22±3 ℃	21.8±0.7 ℃	
Humidity	50±20 %RH	44.9±5.1% RH	
Luminous intensity	150~300 Lux	284 Lux	
Noise	60 db less than	49.2 dB	
Ammonia	15 ppm less than	5 ppm less than	

It is authenticated that there is no change of environment which digresses from the above established value for more than 2 hours during the test period.

Facility management director

Dong-Seok Beck

2013-07-31

Certific	cation of Environment for animal breeding room
Study No.	GT13-00019
Title	Mammalian Erythrocyte Micronucleus Test of MWCNT in ICR mice (Main test)
SPF Room No.	SPF #4 Animal Room
Period of animal Breeding	2014-05-14 ~ 2014-05-22

Breeding environment condition

Section	Range of SOP	Survey value	Remark
Temperature	22±3 ℃	21.5±1.1 ℃	
Humidity	50±20 %RH	53.8±5.1% RH	
Luminous intensity	150~300 Lux	284 Lux	
Noise	60 db less than	49.2 dB	
Ammonia	15 ppm less than	5 ppm less than	


It is authenticated that there is no change of environment which digresses from the above established value for more than 2 hours during the test period.

Facility management director

Dong-Seok Beck

2013-07-31

Annex 8. GLP Certificate



			(뒤 쪽)	-1
< 변	경사항>			
	일자	내용	확인	
200	9.5.20	S명책임자 변경 : 유일 자 (II-Je Yu) 에서 용경 석 (Kyung-Seuk Sorg)으로 변경	E.	
300	(者)	(स्थि) है।	(2)	
	(85)	Test Sage: Partie dermal texacity, Fish: outle texacity (balidation: Since Nov. 16,2009). 2月至 树: 2分至 (Chang-ro Kim) 中心	E	
2010	1.8.2	과표자 변경 : 신청호 (Chang-re Kim) 아시	G L P	
20/0.	8.2	가한 변경 : 한국건당성활환성(PM 파이타왕당보고)'로 변경 기한병 변경 : 한국건당성활환성(PM 파이타왕당보고)'로 변경	G L P	
2011	. 9.9	* 영화점 (Bicconvergence Technology Division, Korea Confirmi 영화업과 변경 : 용 경 설 (Kyung Seuk Song) 에서 이 전 귀 (Jin Kyu Lee) 또한 변경	Ty Laboratories	
<처	분사항>	ा राजा (Uni Ngu Lee) <u>१५५ वर्ष</u>	न्त्र थ	
	일자	내용	확인	
	-			
<참:	고사항>			
	일자	પાક	확인	
	0.12.	전가시주락가 결과, GLP가정을 주수하고 있음	GLP	
	2. 7. 2	(中P Compliance) る水中男子 多次。 Capting 子谷と 文語	확 인 G L P	
207	,,,	(GLT Compliance)	화 인	
	-			

일자 내용 확인 2011. 9. 9 기관업변경: "한국건설생활환경시험연구원 바이오용합단"으로 변경 (Boconvergence Technology Department, Korea Conformity Laboratories) 되고 한경 : 오 대스 (Ticeshik Oh.) 에서 (C. L. P. 수 자신 (Ticeshik Oh.) 에서 함 인 기관정변경: '현지산영광건상-협관위원 파이오등관유건'로 변형 (다. P. P. 기관정변경: '현지산영광건상-협관위원 파이오등관유건'로 변형 (다. P. P.	
2011. 9. 9 기관명변경: "한국건설생활환경시원연구원 바이오용합단"으로 변경 (Boconvergence Technology Department, Korea Conformity Laboratories) 확 인 과도자 변경 : 오 패스 (Taesch)k Oh) 에서 G L P	
(Bicconvergence Technology Department, Korea Conformity Laboratories) 型 인 대도가 변경: 오퍼스 (Taesh)k (Dh) 에서 G L P	
대도자 변경 : 오퍼스 (Taeshik Oh) 예상 GLP	
) (
2012. 7. 2 (Biocontergence Technica as Laboratory Korea Contrainty Talyantia)	
大学 中部: 多世界 子の音がは、大学をなるないが、 to be content to the content to to be refer to the content to the c	P
2012. 7. 2 Test Scape: Daylinia sp: acute toricity, Algae: growth inhibition (since July. 2. 2012)	인
<처분사항>	
일자 내용 확인	
40	
	_
	_
<참고사항>	
일자 내용 확인	

Annex 9. Quality assurance statement-Original

신뢰성보증확인서

시험번호: GT13-00019

시 험 명 : ICR 마우스 골수세포를 이용한 MWCNT의 소핵시험

이 보고서에 기술된 시험을 독립적으로 아래와 같이 시험과정 단계별로 점검하였으며 각 점검결과를 표준작업지원서에 따라 시험책임자와 운영책임자에게 통보 및 보고하였다.

본 시험은 국립환경과학원 고시 제2012-23호(2012년 08월 22일) '화학물칠 유해성시험연구기관의 지정 등에 관한 규정', 동규정 별표5 제4장 제16항 유전독 성시험(골수세포 소핵시험) 및 OECD Guidelines for the Testing of Chemical No. 474에 따라 수행되었으며, 보고서의 작성 및 결과의 기술이 시험의 실시과정에 서 발생한 시험기초자료를 바탕으로 정확히 반영되었음을 확인하였다.

점검내용	실 시 일	시험책임자에게 통보일	운영책임자에게 보고일
시험계획서	2013, 03, 18	2013, 03, 18	2013. 03. 18
동물업수	2013. 05. 14	2013, 05, 14	2013, 05, 14
시험물질 및 대조물질	2013. 04. 01	2013, 04, 01	2013. 04. 01
시험물질조제	2013, 05, 21	2013, 05, 21	2013. 05. 21
동물사육 및 투여	2013. 05: 21	2013. 05. 21	2013. 05. 21
증상관찰 및 측정	2013, 05, 21	2013. 05. 21	2013. 05. 21
부검(골수) 및 검제제작	2013. 05. 22	2013, 05, 22	2013, 05, 22
건경	2013. 05. 22	2013, 05, 22	2013. 05. 22
시험기초자료	2013, 08, 23	2013. 08. 23	2013, 08, 23
최종보고서	2013. 08. 23	2013, 08, 23	2013, 08, 23



한국건설생활환경시험연구원 바이모융합연구소 신뢰성보증책임자 중 경사 (20)

2013년 0년 월구국일

Annex 10. Study personnel-Original

시험관계자서명

