

84-2 Salmonella Mutagenicity Test

Reviewed by: John H.S. Chen
Section I, Toxicology Branch (TS-769C)
Tertiary Reviewer: I. Mauer
Section VI, Toxicology Branch (TS-769C)
Reviewed by Section Head: R.B. Jaeger
Section I, Toxicology Branch (TS-769C)

DATA EVALUATION REPORT

Study Type: Gene Mutation in Bacteria

TOX. CHEM. NO.: 2980

Accession No.: 403888-04

MRID NO.:

Test Material: CGA-154281 Technical (Batch No. HM 4287; 99.3% Purity)

Synonyms:

Study Number (s): 840232

Sponsor: CIBA-GEIGY Corporation Agricultural Division

Testing Facility: CIBA-GEIGY Limited Experimental Pathology Laboratory, Basle, Switzerland

Title of Report: Salmonella/Mammalian-Microsome Mutagenicity Test
(Test material: CGA 154281 technical)

Author(s): E. Deparade and P. Arni

Report Issued: June 14, 1984

Conclusions:

CGA-154281 technical was nonmutagenic to TA98, TA100, TA1535 and TA1537 strains of Salmonella typhimurium either with or without metabolic activation at the concentrations tested.

Concentrations tested: 20, 80, 320, 1280 and 5120 ug/plate

Classification of Data: Unacceptable
(Deficiencies are identified in the detailed review)

Title of Study: Salmonella/Mammalian-Microsome Mutagenicity Test with
CGA-154281 Technical
Giba-Geigy Limited Experimental Pathology Test No. 840232

I. Materials and Methods:

1. Test Materials:

The test compound, CGA-154281 technical (Batch No. HM 4287; 99.3% Purity), was dissolved in acetone. Solutions of daunorubicin-HCl in phosphate buffer, 4-nitroquinoline-N-oxide in phosphate buffer, N-methyl-N'-nitro-N-nitrosoguanidine in phosphate buffer, 9-aminoacridine hydrochloride in DMSO, cyclophosphamide in phosphate buffer and 2-aminoanthracene in DMSO were prepared prior to use and served as positive controls.

2. Bacteria

Four histidine-auxotrophic strains of Salmonella typhimurium (TA98, TA100, TA1535 and TA1537) originally obtained from Dr. Ames were used in this study.

3. In Vitro Metabolic Activation System

The mammalian metabolic activation system consisted of rat liver homogenate from Aroclor 1254-treated rats (Tif:RAIF(SPF)) and a solution of co-factors described by Ames et al. (Mutation Res., 31: 347-364, 1975).

4. Mutagenicity Test

The mutagenicity tests were carried out in accordance with the method described by Ames et al. (1975). The mutagenicity of CGA-154281 technical was evaluated by the Ames test at the concentrations of 20, 80, 320, 1280 and 5120 ug/plate either in the presence or absence of metabolic activation. Mutations were quantified on triplicate plates for each strain by counting the his⁺ revertant colonies after 48 hours of incubation at 37°C on a selective agar plate. Positive control compounds and negative (solvent) control were run concurrently with the test compound.

II. Reported Results:(Tables 1, 2, 3, 4, 5, 6, 7 and 8 attached)

1. Primary Toxicity Test

From the results obtained (Tables 1, 2, 5 and 6), a reduction in the revertant colonies was observed in the CGA-154281-treated cultures of TA1537 strain at the concentration of 5120 ug/plate. Also at the concentration of 5120 ug/plate, the material precipitated in soft agar. Therefore, the highest concentration suitable for the mutagenicity test was found to be 5120 ug/plate.

2

2. Mutagenicity Test

No increase in the number of revertant colonies (less than 2-fold) over concurrent control value was observed for any of tester strains following exposure to the test compound (i.e., 20, 80, 320, 1280 and 5120 ug/plate) either in the presence or absence of metabolic activation.

III. Evaluation and Recommendation:

1. The specific procedures used for confirming the genotypes of TA98, TA100, TA1535, and TA1537 strains of Salmonella typhimurium in accordance with the individual sensitivity test recommended by the Ames test were not given in the report.
2. The spontaneous revertant colonies for each of the four tester strains of Salmonella typhimurium are found within the normal ranges of revertant colonies recommended by the Ames test (Mutation Res., 31: 347-364, 1975).
3. The strain specific control compounds (daunorubicin-HCl; 4-nitroquinoline-N-oxide; N-methyl-N'-nitro-N-nitrosoguanidine; 9-aminoacridine-HCl) and the positive controls (cyclophosphamide; 2-aminocanthracene) to ensure the efficacy of the activation system have given the strong positive responses as expected.
4. Since the cytotoxicity of the test compound against TA1537 strain was observed at 5120 ug/plate in the experiments performed with and without metabolic activation (See results in Tables 1, 2, 5 and 6), the highest dose of CGA-154281 technical (5120 ug/plate) used in this study is considered acceptable. However, the assay should be evaluated using a minimum of five concentrations with adequate intervals between test points (i.e., a narrow range of concentrations recommended). It is, therefore, questionable whether appropriate median doses of the test compound (i.e., 2000, 3000 or 4000 ug/plate) were chosen for this study. The study is judged unacceptable in the present form.

4-(Dichloroacetyl)-3,4-Dihydro-3-Methyl-
2H-1,4-Benzoxazine

Page _____ is not included in this copy.

Pages 4 through 11 are not included.

The material not included contains the following type of information:

- ☐ Identity of product inert ingredients.
 - ☐ Identity of product impurities.
 - ☐ Description of the product manufacturing process.
 - ☐ Description of quality control procedures.
 - ☐ Identity of the source of product ingredients.
 - ☐ Sales or other commercial/financial information.
 - ☐ A draft product label.
 - ☐ The product confidential statement of formula.
 - ☐ Information about a pending registration action.
 - ☒ FIFRA registration data.
 - ☐ The document is a duplicate of page(s) _____.
 - ☐ The document is not responsive to the request.
-

The information not included is generally considered confidential by product registrants. If you have any questions, please contact the individual who prepared the response to your request.
