

National Functional Guidelines Report #03

Lab MITKEM(Mitkem Laboratories) SDG MF2A29 Case 43795 Contract EPW09039 Region 6 DDTID 183441 SOW ISM01.3

Data Review Reports

Blanks

Blanks	ICP_AES
ND03	<p>The following samples have analyte results greater than or equal to MDLs but less than CRQLs. The associated ICB analyte results are greater than or equal to MDLs but less than or equal to CRQLs. Detected analytes are qualified U. Nondetected analytes are not qualified. Sample results are elevated to CRQLs.</p> <p>MF2A53, MF2A54, MF2A61, MF2A62, MF2A63, MF2A33, MF2A34, MF2A48, MF2A50, MF2A52, MF2A55, MF2A56, MF2A57, MF2A57D, MF2A57L, MF2A58, MF2A59, MF2A60, MF2A73, MF2A75, PBS50, MF2A29</p> <p>Sodium MF2A53 , MF2A54 , MF2A61 , MF2A62 , MF2A63 , MF2A33 , MF2A34 , MF2A48 , MF2A50 , MF2A52 , MF2A55 , MF2A56 , MF2A57 , MF2A57D , MF2A57L , MF2A58 , MF2A59 , MF2A60 , MF2A73 , MF2A75 , PBS50</p> <p>Calcium MF2A53 , MF2A54 , MF2A63 , MF2A50 , MF2A52 , MF2A57 , MF2A57L , MF2A60</p> <p>Potassium MF2A53 , MF2A54 , MF2A63 , MF2A29 , MF2A33 , MF2A48 , MF2A50 , MF2A52 , MF2A57 , MF2A57D , MF2A57L , MF2A59 , MF2A60 , PBS50</p>
ND04	<p>The following samples have analyte results greater than or equal to MDLs but less than CRQLs. The associated CCB analyte results are greater than or equal to MDLs but less than or equal to CRQLs. Detected analytes are qualified U. Nondetected analytes are not qualified. Sample results are elevated at CRQLs.</p> <p>MF2A53, MF2A54, MF2A61, MF2A62, MF2A63, MF2A33, MF2A34, MF2A48, MF2A50, MF2A52, MF2A55, MF2A56, MF2A57, MF2A57D, MF2A57L, MF2A58, MF2A59, MF2A60, MF2A73, MF2A75, PBS50, MF2A29</p> <p>Sodium MF2A53 , MF2A54 , MF2A61 , MF2A62 , MF2A63 , MF2A33 , MF2A34 , MF2A48 , MF2A50 , MF2A52 , MF2A55 , MF2A56 , MF2A57 , MF2A57D , MF2A57L , MF2A58 , MF2A59 , MF2A60 , MF2A73 , MF2A75 , PBS50</p> <p>Calcium MF2A53 , MF2A54 , MF2A63 , MF2A50 , MF2A52 , MF2A57 , MF2A57L , MF2A60</p> <p>Potassium MF2A53 , MF2A54 , MF2A63 , MF2A29 , MF2A33 , MF2A48 , MF2A50 , MF2A52 , MF2A57 , MF2A57D , MF2A57L , MF2A59 , MF2A60 , PBS50</p>
ND05	<p>The following samples have analyte results greater than CRQLs. The associated ICB analyte results are greater than or equal to MDLs but less than or equal to CRQLs. Use professional judgment to qualified detected and nondetected analytes.</p> <p>LCS50, MF2A29, MF2A37, MF2A61, MF2A62, MF2A33, MF2A34, MF2A48, MF2A55, MF2A56, MF2A57D, MF2A58, MF2A59, MF2A73, MF2A75</p> <p>Sodium LCS50 , MF2A29 , MF2A37</p> <p>Calcium MF2A61 , MF2A62 , LCS50 , MF2A29 , MF2A33 , MF2A34 , MF2A37 , MF2A48 , MF2A55 , MF2A56 , MF2A57D , MF2A58 , MF2A59 , MF2A73 , MF2A75</p> <p>Potassium MF2A61 , MF2A62 , LCS50 , MF2A34 , MF2A37 , MF2A55 , MF2A56 , MF2A58 , MF2A73 , MF2A75</p>
ND06	<p>The following samples have analyte results greater than CRQLs. The associated CCB analyte results are greater than or equal to MDLs but less than or equal to CRQLs. Use professional judgment to qualified detected and nondetected analytes.</p> <p>LCS50, MF2A37, MF2A61, MF2A62, MF2A29, MF2A33, MF2A34, MF2A48, MF2A55, MF2A56, MF2A57D, MF2A58, MF2A59, MF2A73, MF2A75</p> <p>Sodium LCS50 , MF2A37</p> <p>Calcium MF2A61 , MF2A62 , LCS50 , MF2A29 , MF2A33 , MF2A34 , MF2A37 , MF2A48 , MF2A55 , MF2A56 , MF2A57D , MF2A58 , MF2A59 , MF2A73 , MF2A75</p> <p>Potassium MF2A61 , MF2A62 , LCS50 , MF2A34 , MF2A37 , MF2A55 , MF2A56 , MF2A58 , MF2A73 , MF2A75</p>
NE03	<p>The following samples have analyte results greater than CRQLs but less than 10x preparation blank results. The associated preparation blank analyte results are greater than CRQLs. Detected analytes are qualified J+. Nondetected analytes are not qualified.</p>

National Functional Guidelines Report #03

Lab MITKEM(Mitkem Laboratories) SDG MF2A29 Case 43795 Contract EPW09039 Region 6 DDTID 183441 SOW ISM01.3

Data Review Reports

Blanks

Blanks	ICP_AES
	LCS50
	Iron LCS50
Blanks	ICP_AES
NE04	The following samples have analyte results greater than or equal to MDLs but less than or equal to CRQLs. The associated preparation blank analyte results are greater than or equal to MDLs but less than or equal to CRQLs. Detected analytes are qualified U. Nondetected analytes are not qualified. Sample results are elevated to CRQLs.
	MF2A52, MF2A53, MF2A54, MF2A60, MF2A61, MF2A62, MF2A33, MF2A34, MF2A48, MF2A50, MF2A55, MF2A56, MF2A57, MF2A57D, MF2A58, MF2A59, MF2A63, MF2A73, MF2A75, MF2A29
	Sodium MF2A52 , MF2A53 , MF2A54 , MF2A60 , MF2A61 , MF2A62 , MF2A33 , MF2A34 , MF2A48 , MF2A50 , MF2A55 , MF2A56 , MF2A57 , MF2A57D , MF2A58 , MF2A59 , MF2A63 , MF2A73 , MF2A75
	Potassium MF2A53 , MF2A54 , MF2A60 , MF2A29 , MF2A33 , MF2A48 , MF2A50 , MF2A52 , MF2A57 , MF2A57D , MF2A59 , MF2A63
Blanks	ICP_AES
NE05	The following samples have analyte results greater than CRQLs. The associated preparation blank analyte results are greater than or equal to MDLs but less than or equal to CRQLs. Use professional judgment to qualify detected and nondetected analytes.
	LCS50, MF2A29, MF2A37, MF2A61, MF2A62, MF2A34, MF2A55, MF2A56, MF2A58, MF2A73, MF2A75
	Sodium LCS50 , MF2A29 , MF2A37
	Potassium MF2A61 , MF2A62 , LCS50 , MF2A34 , MF2A37 , MF2A55 , MF2A56 , MF2A58 , MF2A73 , MF2A75

National Functional Guidelines Report #03

Lab MITKEM(Mitkem Laboratories) SDG MF2A29 Case 43795 Contract EPW09039 Region 6 DDTID 183441 SOW ISM01.3

Data Review Reports

Blanks

Blanks	ICP_MS
ND03	<p>The following samples have analyte results greater than or equal to MDLs but less than CRQLs. The associated ICB analyte results are greater than or equal to MDLs but less than or equal to CRQLs. Detected analytes are qualified U. Nondetected analytes are not qualified. Sample results are elevated to CRQLs.</p> <p>MF2A56, MF2A29, MF2A33, MF2A60, MF2A48, MF2A50, MF2A52, MF2A53, MF2A54, MF2A57, MF2A57D, MF2A57L, MF2A59, MF2A63, PBS53, MF2A34, MF2A37, MF2A55, MF2A61, MF2A62, MF2A58, MF2A73, MF2A75</p> <p>Vanadium MF2A56 , MF2A29 , MF2A33 , MF2A60 , MF2A48 , MF2A50 , MF2A52 , MF2A53 , MF2A54 , MF2A57 , MF2A57D , MF2A57L , MF2A59 , MF2A63</p> <p>Cobalt MF2A56 , MF2A29 , MF2A33 , MF2A48 , MF2A50 , MF2A52 , MF2A54 , MF2A57 , MF2A57D , MF2A59</p> <p>Barium PBS53 , MF2A60 , MF2A52 , MF2A53 , MF2A54 , MF2A63 , MF2A57L</p> <p>Beryllium MF2A34 , MF2A37 , MF2A55 , MF2A29 , MF2A61 , MF2A62 , MF2A58 , MF2A73 , MF2A75</p> <p>Thallium MF2A29</p> <p>Cadmium MF2A34 , MF2A37 , MF2A55 , MF2A29 , MF2A61 , MF2A62 , MF2A58 , MF2A73 , MF2A75</p> <p>Silver MF2A75 , MF2A34 , MF2A37 , MF2A55 , MF2A29 , MF2A60 , MF2A48 , MF2A50 , MF2A52 , MF2A54 , MF2A61 , MF2A63</p> <p>Lead MF2A57L</p>
ND04	<p>The following samples have analyte results greater than or equal to MDLs but less than CRQLs. The associated CCB analyte results are greater than or equal to MDLs but less than or equal to CRQLs. Detected analytes are qualified U. Nondetected analytes are not qualified. Sample results are elevated at CRQLs.</p> <p>MF2A56, MF2A29, MF2A33, MF2A60, MF2A48, MF2A50, MF2A52, MF2A53, MF2A54, MF2A57, MF2A57D, MF2A57L, MF2A59, MF2A63, PBS53, MF2A34, MF2A37, MF2A55, MF2A61, MF2A62, MF2A58, MF2A73, MF2A75</p> <p>Vanadium MF2A56 , MF2A29 , MF2A33 , MF2A60 , MF2A48 , MF2A50 , MF2A52 , MF2A53 , MF2A54 , MF2A57 , MF2A57D , MF2A57L , MF2A59 , MF2A63</p> <p>Cobalt MF2A56 , MF2A29 , MF2A33 , MF2A48 , MF2A50 , MF2A52 , MF2A54 , MF2A57 , MF2A57D , MF2A59</p> <p>Barium PBS53 , MF2A60 , MF2A52 , MF2A53 , MF2A54 , MF2A63 , MF2A57L</p> <p>Beryllium MF2A34 , MF2A37 , MF2A55 , MF2A29 , MF2A61 , MF2A62 , MF2A58 , MF2A73 , MF2A75</p> <p>Thallium MF2A29</p> <p>Cadmium MF2A75 , MF2A34 , MF2A37 , MF2A55 , MF2A29 , MF2A61 , MF2A62 , MF2A58 , MF2A73 , MF2A75</p> <p>Copper MF2A33 , MF2A60 , MF2A50 , MF2A52 , MF2A53 , MF2A54 , MF2A62 , MF2A57 , MF2A57D , MF2A59</p> <p>Silver MF2A75 , MF2A34 , MF2A37 , MF2A55 , MF2A29 , MF2A60 , MF2A48 , MF2A50 , MF2A52 , MF2A54 , MF2A61 , MF2A63</p> <p>Lead MF2A57L</p>
ND05	<p>The following samples have analyte results greater than CRQLs. The associated ICB analyte results are greater than or equal to MDLs but less than or equal to CRQLs. Use professional judgment to qualified detected and nondetected analytes.</p> <p>MF2A75, MF2A34, MF2A37, MF2A55, LCS53, MF2A61, MF2A62, MF2A57S, MF2A58, MF2A73, MF2A56, MF2A29, MF2A33, MF2A48, MF2A50, MF2A57, MF2A57D, MF2A59, MF2A57A, MF2A60, MF2A52, MF2A53, MF2A54, MF2A63</p>

National Functional Guidelines Report #03

Lab MITKEM(Mitkem Laboratories) SDG MF2A29 Case 43795 Contract EPW09039 Region 6 DDTID 183441 SOW ISM01.3

Data Review Reports

Blanks

Blanks	ICP_MS
	Vanadium MF2A75 , MF2A34 , MF2A37 , MF2A55 , LCS53 , MF2A61 , MF2A62 , MF2A57S , MF2A58 , MF2A73
	Barium MF2A37 , MF2A55 , MF2A56 , LCS53 , MF2A29 , MF2A33 , MF2A34 , MF2A61 , MF2A48 , MF2A50 , MF2A62 , MF2A57 , MF2A57D , MF2A57S , MF2A58 , MF2A59 , MF2A73 , MF2A75
	Cobalt MF2A75 , MF2A34 , MF2A37 , MF2A55 , LCS53 , MF2A61 , MF2A62 , MF2A57S , MF2A58 , MF2A73
	Beryllium LCS53 , MF2A57S
	Antimony LCS53 , MF2A57A , MF2A57S
	Thallium LCS53 , MF2A57S
	Cadmium LCS53 , MF2A57S
	Silver LCS53 , MF2A57S
	Lead MF2A75 , MF2A34 , MF2A37 , MF2A55 , MF2A56 , LCS53 , MF2A29 , MF2A33 , MF2A60 , MF2A48 , MF2A50 , MF2A52 , MF2A53 , MF2A54 , MF2A61 , MF2A62 , MF2A63 , MF2A57 , MF2A57D , MF2A57S , MF2A58 , MF2A59 , MF2A73
Blanks	ICP_MS
ND06	The following samples have analyte results greater than CRQLs. The associated CCB analyte results are greater than or equal to MDLs but less than or equal to CRQLs. Use professional judgment to qualified detected and nondetected analytes.
	MF2A75, MF2A34, MF2A37, MF2A55, LCS53, MF2A61, MF2A57S, MF2A58, MF2A73, MF2A56, MF2A29, MF2A33, MF2A48, MF2A50, MF2A62, MF2A57, MF2A57D, MF2A57L, MF2A59, MF2A57A, MF2A60, MF2A52, MF2A53, MF2A54, MF2A63
	Vanadium MF2A75 , MF2A34 , MF2A37 , MF2A55 , LCS53 , MF2A61 , MF2A57S , MF2A58 , MF2A73
	Barium MF2A34 , MF2A37 , MF2A55 , MF2A56 , LCS53 , MF2A29 , MF2A33 , MF2A61 , MF2A48 , MF2A50 , MF2A62 , MF2A57 , MF2A57D , MF2A57L , MF2A57S , MF2A58 , MF2A59 , MF2A73 , MF2A75
	Cobalt MF2A75 , MF2A34 , MF2A37 , MF2A55 , LCS53 , MF2A61 , MF2A57S , MF2A58 , MF2A73
	Beryllium LCS53 , MF2A57S
	Antimony LCS53 , MF2A57A , MF2A57S
	Thallium LCS53 , MF2A57S
	Cadmium LCS53 , MF2A57S
	Copper MF2A75 , MF2A48 , MF2A61 , MF2A57S , MF2A58 , MF2A73
	Silver LCS53 , MF2A57S
	Lead MF2A75 , MF2A34 , MF2A37 , MF2A55 , MF2A56 , LCS53 , MF2A29 , MF2A33 , MF2A60 , MF2A48 , MF2A50 , MF2A52 , MF2A53 , MF2A54 , MF2A61 , MF2A62 , MF2A63 , MF2A57 , MF2A57D , MF2A57L , MF2A57S , MF2A58 , MF2A59 , MF2A73
Blanks	ICP_MS
NE04	The following samples have analyte results greater than or equal to MDLs but less than or equal to CRQLs. The associated preparation blank analyte results are greater than or equal to MDLs but less than or equal to CRQLs. Detected analytes are qualified U. Nondetected analytes are not qualified. Sample results are elevated to CRQLs.
	MF2A60, MF2A52, MF2A53, MF2A54, MF2A63

National Functional Guidelines Report #03

Lab MITKEM(Mitkem Laboratories) SDG MF2A29 Case 43795 Contract EPW09039 Region 6 DDTID 183441 SOW ISM01.3

Data Review Reports

Blanks

Blanks	ICP_MS
	Barium MF2A60 , MF2A52 , MF2A53 , MF2A54 , MF2A63
Blanks	ICP_MS
NE05	The following samples have analyte results greater than CRQLs. The associated preparation blank analyte results are greater than or equal to MDLs but less than or equal to CRQLs. Use professional judgment to qualify detected and nondetected analytes.
	MF2A75, MF2A34, MF2A55, MF2A56, LCS53, MF2A29, MF2A33, MF2A59, MF2A37, MF2A48, MF2A50, MF2A61, MF2A62, MF2A57, MF2A57D, MF2A57S, MF2A58, MF2A73
	Barium MF2A75 , MF2A34 , MF2A55 , MF2A56 , LCS53 , MF2A29 , MF2A33 , MF2A59 , MF2A37 , MF2A48 , MF2A50 , MF2A61 , MF2A62 , MF2A57 , MF2A57D , MF2A57S , MF2A58 , MF2A73

National Functional Guidelines Report #03

Lab MITKEM(Mitkem Laboratories) **SDG** MF2A29 **Case** 43795 **Contract** EPW09039 **Region** 6 **DDTID** 183441 **SOW** ISM01.3

Data Review Reports

Detection Limit

Detection Limit	Hg
NDL1	The following samples have results greater than or equal to MDLs but less than CRQLs. Detected analytes are qualified J.
	MF2A29, MF2A48, MF2A50, MF2A52, MF2A54, MF2A55, MF2A56, MF2A57, MF2A57D, MF2A58, MF2A60, MF2A61, MF2A62, MF2A75
	Mercury MF2A29 , MF2A48 , MF2A50 , MF2A52 , MF2A54 , MF2A55 , MF2A56 , MF2A57 , MF2A57D , MF2A58 , MF2A60 , MF2A61 , MF2A62 , MF2A75

National Functional Guidelines Report #03

Lab MITKEM(Mitkem Laboratories) SDG MF2A29 Case 43795 Contract EPW09039 Region 6 DDTID 183441 SOW ISM01.3

Data Review Reports

Detection Limit

Detection Limit	ICP_AES
NDL1	The following samples have results greater than or equal to MDLs but less than CRQLs. Detected analytes are qualified J.
	MF2A52, MF2A53, MF2A60, MF2A61, MF2A62, MF2A33, MF2A34, MF2A48, MF2A50, MF2A54, MF2A55, MF2A56, MF2A57, MF2A57D, MF2A57L, MF2A58, MF2A59, MF2A63, MF2A73, MF2A75, PBS50, MF2A29
	Sodium MF2A52 , MF2A53 , MF2A60 , MF2A61 , MF2A62 , MF2A33 , MF2A34 , MF2A48 , MF2A50 , MF2A54 , MF2A55 , MF2A56 , MF2A57 , MF2A57D , MF2A57L , MF2A58 , MF2A59 , MF2A63 , MF2A73 , MF2A75 , PBS50
	Calcium MF2A53 , MF2A54 , MF2A50 , MF2A52 , MF2A57 , MF2A57L , MF2A60 , MF2A63
	Potassium MF2A52 , MF2A53 , MF2A54 , MF2A60 , MF2A29 , MF2A33 , MF2A48 , MF2A50 , MF2A57 , MF2A57D , MF2A57L , MF2A59 , MF2A63 , PBS50
	Magnesium MF2A54 , MF2A33 , MF2A48 , MF2A50 , MF2A52 , MF2A57 , MF2A57D , MF2A59 , MF2A60 , MF2A63
	Manganese MF2A57L

National Functional Guidelines Report #03

Lab MITKEM(Mitkem Laboratories) SDG MF2A29 Case 43795 Contract EPW09039 Region 6 DDTID 183441 SOW ISM01.3

Data Review Reports

Detection Limit

Detection Limit	ICP_MS
NDL1	The following samples have results greater than or equal to MDLs but less than CRQLs. Detected analytes are qualified J.
	MF2A54, MF2A29, MF2A33, MF2A59, MF2A60, MF2A48, MF2A50, MF2A52, MF2A53, MF2A56, MF2A57, MF2A57D, MF2A57L, MF2A63, MF2A62, MF2A73, PBS53, MF2A75, MF2A34, MF2A55, MF2A37, MF2A61, MF2A58
	Vanadium MF2A54 , MF2A29 , MF2A33 , MF2A59 , MF2A60 , MF2A48 , MF2A50 , MF2A52 , MF2A53 , MF2A56 , MF2A57 , MF2A57D , MF2A57L , MF2A63
	Arsenic MF2A56 , MF2A29 , MF2A33 , MF2A48 , MF2A62 , MF2A57 , MF2A59 , MF2A73
	Chromium MF2A56 , MF2A33 , MF2A59 , MF2A60 , MF2A48 , MF2A50 , MF2A52 , MF2A53 , MF2A54 , MF2A57 , MF2A57D , MF2A63
	Cobalt MF2A29 , MF2A33 , MF2A59 , MF2A48 , MF2A50 , MF2A52 , MF2A54 , MF2A56 , MF2A57 , MF2A57D
	Barium PBS53 , MF2A60 , MF2A52 , MF2A53 , MF2A54 , MF2A57L , MF2A63
	Zinc MF2A33 , MF2A53 , MF2A57 , MF2A57D
	Nickel MF2A54 , MF2A33 , MF2A59 , MF2A60 , MF2A48 , MF2A50 , MF2A52 , MF2A53 , MF2A56 , MF2A57 , MF2A57D , MF2A63
	Beryllium MF2A75 , MF2A34 , MF2A55 , MF2A29 , MF2A37 , MF2A61 , MF2A62 , MF2A58 , MF2A73
	Thallium MF2A29
	Cadmium MF2A75 , MF2A34 , MF2A55 , MF2A29 , MF2A37 , MF2A61 , MF2A62 , MF2A58 , MF2A73
	Copper MF2A33 , MF2A59 , MF2A60 , MF2A50 , MF2A52 , MF2A53 , MF2A54 , MF2A62 , MF2A57 , MF2A57D
	Silver MF2A75 , MF2A34 , MF2A54 , MF2A55 , MF2A29 , MF2A60 , MF2A37 , MF2A48 , MF2A50 , MF2A52 , MF2A61 , MF2A63
	Lead MF2A57L

National Functional Guidelines Report #03

Lab MITKEM(Mitkem Laboratories) SDG MF2A29 Case 43795 Contract EPW09039 Region 6 DDTID 183441 SOW ISM01.3

Data Review Reports

Duplicates

Duplicates	ICP_AES
NI03	The following Duplicate and original sample results are greater than 5xCRQL and RPD is greater than 20. The original sample results are greater than or equal to MDLs. Detected analytes are qualified J. Nondetected analytes are qualified UJ.
	MF2A29, MF2A33, MF2A34, MF2A37, MF2A48, MF2A50, MF2A52, MF2A53, MF2A54, MF2A55, MF2A56, MF2A57, MF2A58, MF2A59, MF2A60, MF2A61, MF2A62, MF2A63, MF2A73, MF2A75
	Iron MF2A57D

National Functional Guidelines Report #03

Lab MITKEM(Mitkem Laboratories) SDG MF2A29 Case 43795 Contract EPW09039 Region 6 DDTID 183441 SOW ISM01.3

Data Review Reports

Holding Times/Preservation

Holding Times/Preservation	Hg
NHT01	The following preserved samples are improperly maintained at temperatures outside the range of 4+/-2 C. Detected analytes with results greater than or equal to MDLs are qualified J-. Use professional judgment to qualify the nondetected analytes.
	MF2A48, MF2A50, MF2A58, MF2A59, MF2A60

National Functional Guidelines Report #03

Lab MITKEM(Mitkem Laboratories) SDG MF2A29 Case 43795 Contract EPW09039 Region 6 DDTID 183441 SOW ISM01.3

Data Review Reports

Holding Times/Preservation

Holding Times/Preservation	ICP_AES
NHT01	The following preserved samples are improperly maintained at temperatures outside the range of 4+/-2 C. Detected analytes with results greater than or equal to MDLs are qualified J-. Use professional judgment to qualify the nondetected analytes.
	MF2A48, MF2A50, MF2A58, MF2A59, MF2A60

National Functional Guidelines Report #03

Lab MITKEM(Mitkem Laboratories) SDG MF2A29 Case 43795 Contract EPW09039 Region 6 DDTID 183441 SOW ISM01.3

Data Review Reports

Holding Times/Preservation

Holding Times/Preservation	ICP_MS
NHT01	The following preserved samples are improperly maintained at temperatures outside the range of 4+/-2 C. Detected analytes with results greater than or equal to MDLs are qualified J-. Use professional judgment to qualify the nondetected analytes.
	MF2A60, MF2A48, MF2A50, MF2A58, MF2A59

National Functional Guidelines Report #03

Lab MITKEM(Mitkem Laboratories) SDG MF2A29 Case 43795 Contract EPW09039 Region 6 DDTID 183441 SOW ISM01.3

Data Review Reports

Matrix Spikes

Matrix Spikes	ICP_MS
NG11	The following Matrix Spike samples have percent recoveries in the range of 30-74% and post-digestion spike samples have percent recoveries greater than or equal to 75%. Detected analytes with results greater than or equal to MDLs are qualified J. Nondetected analytes are qualified UJ. MF2A29, MF2A33, MF2A34, MF2A37, MF2A48, MF2A50, MF2A52, MF2A53, MF2A54, MF2A55, MF2A56, MF2A57, MF2A58, MF2A59, MF2A60, MF2A61, MF2A62, MF2A63, MF2A73, MF2A75
	Antimony MF2A57S