2-45 RE . 1. 79



E. I. DU PONT DE NEMOURS & COMPANY INCOMPORATED

WILMINGTON, DELAWARE 19898

POLYMER PRODUCTS DEPARTMENT EXPERIMENTAL STATION

PERSONAL AND CONFIDENTIAL

A. J. Dahl cc: B. W. Karrh L. J. Papa Pral File I.C.

- 269 Complainant's

Exhibit No.

- 353

- N11400

July 16, 1981

DR. D. O. INSLEE SPRUANCE TEXTILE FIBERS RICHMOND, VA.

> ANALYSIS OF BLOOD SAMPLES FOR PERFLUOROOCTANOATE (Job No. 811-601; PRAL Nos. 81-2244-2278; Notebook Nos. E22514, E26238)

As requested in your letter of 5/8/81 to L. J. Papa, the 35 blood samples submitted then have been analyzed for perfluorooctanoate (C_8). Results and sample identification are given in the attached table.

As noted there, the analyses were done using a gas chromatographic method specific for C8 (Lab Method Number ES-567) but results have been reported as ppm F for comparison with total organic fluorine analyses. Precision is $\pm~10\%$ relative standard deviation over most of the concentration range, somewhat less at the lowest values. The lower limit for quantitation is 0.007 ppm F (0.01 ppm perfluorooctanoic acid), with a detection limit of \sim 0.004 ppm which can be distinguished from the reagent background but not well quantitated.

Please contact me (772-4440) or L. J. Papa (772-2745) if you have any questions regarding the analyses. General questions on blood sampling can be directed to J. W. Raines or L. F. Percival.

> 5.5J4Mm2 S. S. Stafford

Attachment jah

Key Words:

Perfluorooctanoic Acid Perfluorocctanoate Blood Aralysis GC

There's a world of things we're doing something about

TABLE I

CONCENTRATION OF PERFLUOROOCTANOATE IN BLOOD (a)

| Sample | | | | GC Analysis | | (b) |
|----------|--------------|-----------|------------|-------------------|--------------|-------|
| PRAL No. | Date Sampled | P.R.No. | Name | Date Analyzed | [Cg], Lg F/g | blood |
| 81-2244 | 5/6/81 | 898 | | 6/11/81 | 0.036 | |
| 81-2245 | 4/27/81 | 1084 | | 6/11/81 | 0.022 | |
| 81-2246 | 4/27/81 | 992 | | 6/11/81 | 0.022 | |
| 81-2247 | 5/6/81 | 357 | | 6/12/81 | 0.030 | |
| 81-2248 | 4/28/81 | 983 | | 6/12/81 | 0.027 | |
| 81-2249 | 5/11/81 | 923 | | 6/12/81 | 0.033 | |
| 81-2250 | 4/28/81 | 1605 | | 6/12/81 | 0.041 | |
| 81-2251 | 5/5/81 | 1544 | | 6/15/81 | 0.012 | |
| 81-2252 | 4/28/81 | 1212 | | 6/15/81 | 0.045 | |
| 81-2253 | 4/27/81 | 946 | | 6/15/81 | 0.086 | |
| 81-2254 | 4/28/81 | 1538 | | 6/15/81 | 0.056 | |
| 81-2255 | 5/7/81 | 3060 | | 6/15/81 | n.d. | |
| 81-2256 | 4/27/81 | 2093 | | 6/15/81 | 0.018 | |
| 81-2257 | 5/6/81 | 2457 | | 6/15/81 | <.007 | |
| 81-2258 | 5/7/81 | 1908 | | 6/15/81 | 0.009 | |
| 81-2259 | 5/7/81 | 2542 | | 6/15/81 | n.d. | |
| 81-2260 | 5/14/81 | 2151 | | 6/16/81 | 0.016 | |
| 81-2261 | 5/14/81 | 2205 | | 6/16/81 | 0.011 | * * |
| 81-2262 | 5/15/81 | 651 | | 6/16/81 | n.d. | |
| 81-2263 | 5/7/81 | 3184 | | 6/18/81 | n.d. | |
| 81-2264 | 5/18/81 | 1190 | | 6/17/81 | 0.082 | |
| 81-2265 | 5/15/81 | sal. | | 6/11/81 | n.d. | |
| 81-2266 | 5/15/81 | 1949 | | 6/11/81 | n.d. | |
| 81-2267 | 5/14/81 | 1856 | 有多数 | 6/11/81 | n.d. | |
| 81-2268 | 5/6/81 | 2771 | | 6/12/81 | n.d. | |
| 81-2269 | 5/14/81 | sa1 | | 6/12/81 | n.d. | |
| 81-2270 | 5/7/81 | 1466 | | 6/12/81 | n.d. | |
| 81-2271 | 5/7/81 | sal | | 6/12/81 | n.d. | |
| 81-2272 | 5/7/81 | sal | | 6/15/81 | n.d. | |
| 81-2273 | 5/14/81 | pensioned | | 6/15/81 | n.d. | |
| 81-2274 | 5/7/81 | sal | | 6/15/81 | n.d. | |

TABLE I

CONCENTRATION OF PERFLUOROCCTANOATE IN BLOOD (a)

| Sample | a * *. | | GC Analysis | (b) |
|----------|--------------|--------------|----------------------|--------------------|
| PRAL No. | Date Sampled | P.R.No. Name | Date Analyzed | [Cg], ug F/g blood |
| 81-2275 | 5/7/81 | pensioned | 6/15/81 | n.d. |
| 81-2276 | 5/7/81 | sal | 6/16/81 ¿ 6/18/81 | 0.027 |
| 81-2277 | 5/7/81 | 2898 | 6/16/81 & 6/18/81 | 0.010 |
| 81-2278 | 5/8/81 | pensioned | 6/17/81 | n.d. |

- (a) Analysis as described in Lab Method ES-567 ("Determination of Perfluorooctanoic Acid in Blood, Gas Chromatographic Method", S. Stafford, 4/3/81), using the packed column GC analysis with perfluoro-n-octanoic acid as calibration standard.
- (b) Although the analysis is specifically for perfluorooctanoate (acid or salts), concentrations are given in ppm fluorine for comparison with the results of total organic fluorine analyses. (ppm F = 0.688 x ppm perfluorooctanoic acid) Estimated uncertainty is \pm 10% relative standard deviation. The lower limit for quantitation is 0.007 µgF/g. The detection limit is \sim 0.004 µgF/g, but concentrations in that range cannot be well quantitated and are reported as < 0.007. None detected (n.d.) is reported for samples with [Cg] < 0.004 ppm, which cannot be distinguished from reagent background.