



CONFIDENTIAL

CC: J. W. Karrh - N-11400
 L. J. Papa - 269
 R. W. Ubersax - 353
 PRAL File - 353
 I.C. - 323

5926

E. I. DU PONT DE NEMOURS & COMPANY
 INCORPORATED

WILMINGTON, DELAWARE 19898

POLYMER PRODUCTS DEPARTMENT
 EXPERIMENTAL STATION

PPD

NOV 4 1982

EXPR. STATION
 INFO CENTER

November 1, 1982

Complainant's
 Exhibit No. 90

PERSONAL AND CONFIDENTIAL

TO: J. D. DE LAURENTIS, M.D. - CR&DD, BUILDING 301

FROM: S. R. LAAS

*S.R. Laas*ANALYSIS OF BLOOD SAMPLES FOR PERFLUOROOCTANOATE

(Job No. 820-969, PRAL Nos. 82-4096, 4099, 4100, 5166, 5168, 5172,
 5176-77, 5262, 5269, 5573-75, 5578)

Notebook Nos. E29297, E29298

The 14 blood samples, previously reported in letters to you dated 10/19/82 and 10/29/82 as having perfluorooctanoate (C8) concentrations greater than the blank, have been analyzed for C8 by the usual gas chromatographic method ES-567. Results and sample identifications are given in the attached table.

The samples are all very low; indeed 3 of the samples are indistinguishable from the blank and 5 are below the quantitation limit. The estimated uncertainty for these analyses is higher than usual due to the very low levels of C8 and the difficulty in integrating the small C8 peak which elutes between the 05 and 06 peaks.

If you have any questions regarding these analyses, don't hesitate to call me.

Key Words

Perfluorooctanoate
 GC
 Blood Analysis

JRW000010

SRL:kab
 Attach.

CONFIDENTIAL

TABLE I

CONCENTRATION OF PERFLUOROOCTANOATE IN BLOOD (a)

Sample PRAL No.	Date Received	Name	GC Analysis (b)	
			Date Analyzed	[C8], $\mu\text{g F/g blood}$
82-4096	7/21/82		10/25/82	0.009
82-4099	7/21/82		10/25/82	n.d.
82-4100	7/21/82		10/25/82	<0.007
82-5166	9/17/82		10/22/82	<0.007
82-5168	9/17/82		10/22/82	0.010
82-5172	9/17/82		10/22/82	0.007
82-5176	9/17/82		10/22/82	<0.007
82-5177	9/17/82		10/25/82	0.007
82-5262	9/24/82		10/25/82	<0.007
82-5269	9/24/82		10/25/82	0.009
82-5573	10/8/82		10/26/82	0.007
82-5574	10/8/82		10/26/82	<0.007
82-5575	10/8/82		10/26/82	n.d.
82-5578	10/8/82		10/26/82	n.d.

(a) Analysis as described in Lab Method ES-567 ("Determination of Perfluorooctanoic Acid in Blood, Gas Chromatographic Method", S. Stafford, 4/3/81), using the packed column GC analysis with perfluoro-n-octanoic acid as calibration standard.

(b) Although the analysis is specifically for perfluorooctanoate (acid or salts), concentrations are given in ppm fluorine for comparison with the results of total organic fluorine analyses. ($\text{ppm F} = 0.688 \times \text{ppm perfluorooctanoic acid}$). Estimated uncertainty is $\pm 20\%$ relative standard deviation. The lower limit for quantitation is $0.007 \mu\text{g F/g}$. The detection limit is $\sim 0.004 \mu\text{g F/g}$, but concentrations in that range cannot be well quantitated and are reported as <0.007 . None detected (n.d) is reported for samples with $[\text{C8}] \lesssim 0.004 \text{ ppm}$ which cannot be distinguished from reagent background.

JR W000011