

CONFIDENTIAL

CC: S. C. Croft

I.C.

PRAL File

256 B. W. Culpepper - N-11400

L. J. Papa - 256

256 323

> Complainant's Exhibit No.

E. I. DU PONT DE NEMOURS & COMPANY INCORPORATED

WILMINGTON, DELAWARE 19898

POLYMER PRODUCTS DEPARTMENT EXPERIMENTAL STATION

PPD

September 21, 1983

PERSONAL CONFIDENTIAL

SEP 2 2 1983

EXPR. STATION INFO CENTER

DR. O. P. TEDESCO ENGINEERING DEPARTMENT 11E16 LOUVIERS BUILDING

> ANALYSIS OF BLOOD SAMPLES FOR PERFLUOROOCTANOATE (Job No. 831-1046, PRAL Nos. 83-4848-4849; Notebook Nos. E27433, E29324, E29325)

As requested in your letter of 9/14/83, the two blood samples submitted from Tralee Park (F&FP) employees have been analyzed for perfluorooctanoate (Cg) by the usual gas chromatographic method ES-567. Results and sample identification are given in the attached table.

BETTER THINGS FOR DETTER

Shavon Laas

Sharron Laas

SL:kab Attach.

Key Words

Perfluorooctanoate Blood Analysis

Sample PRAL No.	Date Sampled	P.R.No.	Name	GC Analysis Date Analyzed	[CR], ug F/g blood
83-4848	<u>.</u>	. <u>-</u>		9/16/83	0.17
83-4849		= , ,		9/16/83	0.15

- (a) Analysis as described in Lab Method ES-567 ("Determination of Perfluorocctanoic Acid in Blood, Gas Chromatographic Method", S. Stafford, 4/3/81), using the packed column GC analysis with perfluoro-n-octanoic acid as calibration standard.
- (b) Although the analysis is specifically for perfluorooctanoate (acid or salts), concentrations are given in ppm fluorine for comparison with the results of total organic fluorine analyses. (ppm F = 0.688 x ppm perfluorooctanoic acid) Estimated uncertainty is ± 10% relative standard deviation. The lower limit for quantitation is 0.007 µgF/g. The detection limit is ~ 0.004 µgF/g, but concentrations in that range cannot be well quantitated and are reported as < 0.007. None detected (n.d.) is reported for samples with [C₈] < 0.004 ppm. which cannot be distinguished from reagent background.