



UNITED STATES ENVIRONMENTAL PROTECTION AGENCY
REGION III
Environmental Sciences Center
701 Mapes Road
Fort Meade, Maryland 20755-5350

DATE : February 8, 2012

SUBJECT: Region III Data QA Review

FROM: Colleen Walling *Colleen Walling*
Region III ESAT RPO (3EA20)

TO: Rich Fetzer
Remedial Project Manager (3HS31)

Attached is the inorganic data validation report for the Dimock Residential Groundwater site (Case #: 180-3644-1) completed by the Region III Environmental Services Assistance Team (ESAT) contractor under the direction of Region III EAID.

If you have any questions regarding this review, please call me at (410) 305-2763.

Attachment

TO: #0037 TDF: #02025A

cc: Gene Nance (Techlaw)
Suddha Graves (Techlaw)

OFFICE OF ANALYTICAL SERVICES AND QUALITY ASSURANCE

Printed on 100% recycled/recyclable paper with 100% post-consumer fiber and process chlorine free.
Customer Service Hotline: 1-800-438-2474

Lockheed Martin IS&GS – Civil
Energy & Environment
ESAT Region 3
US EPA Environmental Science Center
701 Mapes Road Ft. Meade, MD 20755-5350

Ex. 4 - CBI

DATE: February 7, 2012

SUBJECT: Level IM2 Inorganic Data Validation for Project 180-3644-1
SDG: H-1
Site: Dimock

FROM: Ex. 4 - CBI
Senior Data Reviewer

Through: Ex. 4 - CBI
Senior Data Review Chemist

TO: Colleen Walling
ESAT Region 3 Project Officer

OVERVIEW

Project 180-3644-1, Sample Delivery Group (SDG) H-1, from the Dimock site consisted of eleven (11) aqueous samples analyzed for total metals by ICP-MS and for mercury (Hg) by cold vapor technique and eleven (11) filtrate samples analyzed for dissolved metals by ICP-MS and for dissolved mercury (Hg) by cold vapor technique. All samples were submitted to Test America Laboratories Pittsburg (TALP). Samples were analyzed for metals and mercury (Hg) in accordance with EPA Methods 6020 and 7470, respectively.

SUMMARY

Data were validated according to the Region 3 Modifications to the National Functional Guidelines for Inorganic Data Review, level IM2 and is assigned the Superfund Data Validation Label S4VM (Stage_4_Validation_Manual). Areas of concern with respect to data usability are listed below.

Data in this Case have been impacted by outliers present in the laboratory blanks as well as the ICP serial dilution analyses. Details for these outliers are discussed under "Minor Problems", specific samples affected are outlined in "Table 1A" and qualified analytical results for all samples are summarized on the Data Summary Forms (DSFs).

MINOR PROBLEMS

Method (MB) and Continuing Calibration (CCB) Blanks had reported results greater than the Method Detection Limits (MDLs) as listed below. Positive results reported for these analytes, which were less than five times ($< 5X$) the blank concentration, may be biased high and have been qualified "B" on the DSFs.

Total Metals:

| <u>Blank</u> | <u>Analytes</u> |
|--------------|---|
| CCB | manganese (Mn), molybdenum (Mo), antimony (Sb) |
| MB | aluminum (Al), arsenic (As), boron (B), chromium (Cr), copper (Cu), Iron (Fe), potassium (K), nickel (Ni), selenium (Se), vanadium (V) |

Dissolved Metals:

| <u>Blank</u> | <u>Analytes</u> |
|--------------|---|
| CCB | Al, cobalt (Co), Mn, Sb, thallium (Tl) |
| MB | silver (Ag), barium (Ba), B, beryllium (Be), Cr, Cu, Fe, lead (Pb), Mo, V, zinc (Zn) |

In the total metals analyses, Percent Differences (%Ds) for the ICP serial dilution analysis were outside the control limit ($>10\%$) for B, magnesium (Mg) and Mn. Positive results for these analytes in affected samples in this matrix are estimated and have been qualified "J" on the DSFs unless superseded by "B".

NOTES

Positive results which are less than the Reporting Limits (RLs) but greater than MDLs have been qualified "J" on the DSFs unless superseded by "B".

Matrix spike, matrix spike duplicate and ICP serial dilution were performed for both total and dissolved samples. Data were qualified accordingly.

In the total metal analyses, several analyte concentrations reported on Form 8 - IN (ICP Serial Dilution) did not match the initial sample results as reported on the Form Is. This from was corrected and the Percent Differences (%Ds) were recalculated by the reviewer. Data were qualified based on the corrected %Ds.

The requirement for ICP-MS and cold vapor MDLs to be performed within one (1) year of sample analysis was not met by the laboratory. No action was taken by the reviewer based on these findings.

ATTACHMENTS

INFORMATION REGARDING REPORT CONTENT

| | |
|-------------------|--|
| TABLES 1A | SUMMARY OF QUALIFIERS ON DATA SUMMARY FORMS AFTER DATA VALIDATION |
| TABLE 1B | CODES USED IN COMMENTS COLUMN OF TABLES 1A |
| APPENDIX A | GLOSSARY OF DATA QUALIFIER CODES |
| APPENDIX B | DATA SUMMARY FORM(S) |
| APPENDIX C | CHAIN OF CUSTODY RECORD(S) |
| APPENDIX D | LABORATORY CASE NARRATIVE(S) |

DCN: Project 180-3644-11M2

TABLE 1A
SUMMARY OF QUALIFIERS ON DATA SUMMARY
FORM AFTER DATA VALIDATION

PROJECT: 180-3644-1

SDG#: H-1

TOTAL METALS

| <u>ANALYTE</u> | <u>SAMPLES AFFECTED</u> | <u>POSITIVE VALUES</u> | <u>NON- DETECTED VALUES</u> | <u>BIAS</u> | <u>COMMENTS*</u> |
|----------------|--|----------------------------|-------------------------------------|-------------|----------------------------------|
| Al | H-1, FH-1, FPT-1, KDE-1, D-1 | B | | High | MB (2.95 J µg/L) |
| As | H-1, FH-1, FPT-1, R-1, R-2, RD-1, RU-1, KDE-1, D-1 | B | | High | MB (0.424 J µg/L) |
| B | All Samples Except D-1 | J | | | ISD (42%) |
| | D-1 | B | | | MB (1.91 J µg/L) ISD (42%) |
| Cr | All Samples | B | | | MB (1.64 J µg/L) |
| Cu | H-1, FH-1, R-1, R-2, RD-1, RU-1, N-1, D-1 | B | | High | MB (0.720 J µg/L) |
| Fe | FPT-1, RD-1, RU-1, KDE-1, N-1, D-1 | B | | High | MB (18.2 J µg/L) |
| K | N-1 | B | | | MB (41.2 J µg/L) |
| Mg | All Samples | J | | | ISD (22%) |
| Mn | All Samples Except N-1, D-1 | J | | | ISD (33%) |
| | N-1, D-1 | B | | | CCB (0.0590 J µg/L) ISD (33%) |

* See explanation of Comments on Table 1B.

TABLE 1A
SUMMARY OF QUALIFIERS ON DATA SUMMARY
FORM AFTER DATA VALIDATION

PROJECT: 180-3644-1

SDG#: H-1

TOTAL METALS

| <u>ANALYTE</u> | <u>SAMPLES AFFECTED</u> | <u>POSITIVE VALUES</u> | <u>NON- DETECTED VALUES</u> | <u>BIAS</u> | <u>COMMENTS*</u> |
|----------------|--|----------------------------|-------------------------------------|-------------|--------------------|
| Ni | H-1, FH-1, FPT-1, R-1, R-2, RD-1, RU-1 | B | | High | MB (0.532 J µg/L) |
| Mo | H-1, FPT-1, S-1, R-1, R-2, RD-1, KDE-1 | B | | High | CCB (0.503 J µg/L) |
| Sb | All Samples | B | | High | CCB (0.704 J µg/L) |
| Se | All Samples | B | | High | MB (1.72 J µg/L) |
| V | FPT-1, R-1, RD-1 D-1 | B | | High | MB (0.492 J µg/L) |

* See explanation of Comments on Table 1B

TABLE 1A
SUMMARY OF QUALIFIERS ON DATA SUMMARY
FORM AFTER DATA VALIDATION

PROJECT: 180-3644-1
SDG#: H-1
DISSOLVED METALS

| <u>ANALYTE</u> | <u>SAMPLES</u> <u>AFFECTED</u> | <u>POSITIVE</u> <u>VALUES</u> | <u>NON-</u> <u>DETECTED</u> <u>VALUES</u> | <u>BIAS</u> | <u>COMMENTS*</u> |
|----------------|---|----------------------------------|---|-------------|---------------------|
| Ag | H-1, FH-1 | B | | High | MB (0.0820 J µg/L) |
| Al | H-1, FH-1, S-1, R-2, RD-1, RU-1, KDE-1 | B | | High | CCB (5.98 J µg/L) |
| Ba | N-1, D-1 | B | | High | MB (0.752 J µg/L) |
| B | FH-1, RD-1, RU-1, N-1, D-1 | B | | | MB (3.03 J µg/L) |
| Be | S-1 | B | | High | MB (0.0370 J µg/L) |
| Co | FH-1, FTP-1, S-1, R-1, R-2, RD-1 | B | | | CCB (0.0430 J µg/L) |
| Cr | All Samples | B | | | MB (1.45 J µg/L) |
| Cu | All Samples Except, FPT-1 | B | | High | MB (0.644 J µg/L) |
| Fe | H-1, FH-1, FPT-1, R-2 | B | | High | MB (19.3 J µg/L) |
| Mn | H-1, FPT-1, S-1, RD-1, RU-1, N-1, D-1 | B | | | CCB (0.409 J µg/L) |
| Pb | All Samples | B | | | MB (0.0690 J µg/L) |
| Mo | H-1, FH-1, FPT-1, S-1, R-1, R-2 | B | | | MB (0.907 J µg/L) |

* See explanation of Comments on Table 1B.

TABLE 1A
SUMMARY OF QUALIFIERS ON DATA SUMMARY
FORM AFTER DATA VALIDATION

PROJECT: 180-3644-1
SDG#: H-1
DISSOLVED METALS

| <u>ANALYTE</u> | <u>SAMPLES AFFECTED</u> | <u>POSITIVE VALUES</u> | <u>NON- DETECTED VALUES</u> | <u>BIAS</u> | <u>COMMENTS*</u> |
|----------------|--|----------------------------|-------------------------------------|-------------|---------------------------------|
| Sb | All Samples | B | | High | CCB (1.18 J $\mu\text{g/L}$) |
| Tl | All Samples Except KDE-1, D-1 | B | | High | CCB (0.0500 J $\mu\text{g/L}$) |
| V | FPT-1, S-1, RU-1 | B | | High | MB (0.119 J $\mu\text{g/L}$) |
| Zn | H-1, FH-1, R-1, RD-1, RU-1, KDE-1 N-1, D-1 | B | | | MB (1.74 J $\mu\text{g/L}$) |

* See explanation of Comments on Table 1B.

TABLE 1B
CODES USED IN COMMENTS COLUMN

- MB = Method blanks (MBs) had results greater than the MDLs (results are in parenthesis). Positive results less than five times (<5X) the value of the blank may be biased high.
- ISD = Percent Differences (%Ds) for the ICP serial dilution analyses exceeded the control limit (10%) [% Ds are in parenthesis]. Reported results are estimated.
- CCB = Continuing calibration blanks (CCBs) had results greater than the MDLs (results are in parenthesis). Positive results less than five times (<5X) the value of the blank may be biased high.

APPENDIX A
Glossary of Qualifier Codes

GLOSSARY OF DATA QUALIFIER CODES (INORGANIC)

CODES RELATED TO IDENTIFICATION

(confidence concerning presence or absence of analytes):

U = Not detected. The associated number indicates approximate sample concentration necessary to be detected.

(NO CODE) = Confirmed identification.

B = Not detected substantially above the level reported in laboratory or field blanks.

R = Unusable result. Analyte may or may not be present in the sample. Supporting data necessary to confirm result.

N = Tentative identification. Consider present.
Special methods may be needed to confirm its presence or absence in future sampling efforts.

CODES RELATED TO QUANTITATION

(can be used for both positive results and sample quantitation limits):

J = Analyte Present. Reported value may not be accurate or precise.

K = Analyte present. Reported value may be biased high. Actual value is expected to be lower.

L = Analyte present. Reported value may be biased low.
Actual value is expected to be higher.

UJ = Not detected, quantitation limit may be inaccurate or imprecise.

UL = Not detected, quantitation limit is probably higher.

OTHER CODES

Q = No analytical result.

APPENDIX B
Data Summary Forms

DATA SUMMARY FORM: INORGANIC

Page 1 of 6

Project #: 180-3644-1

SDG: H-1

Number of Soil Samples : 0

Site : DIMOCK

Number of Water Samples : 11

Lab. : TALP

TOTAL METALS

| | | | | | | | | | | | |
|-------------------|------|------------|------|------------|------|------------|------|------------|------|------------|------|
| Sample Number: | | H-1 | | FH-1 | | FPT-1 | | S-1 | | R-1 | |
| Station Location: | | H-1 | | FH-1 | | FPT-1 | | S-1 | | R-1 | |
| Laboratory ID: | | 180-3644-1 | | 180-3644-2 | | 180-3644-3 | | 180-3644-4 | | 180-3644-5 | |
| Field QC: | | | | | | | | | | | |
| Matrix : | | Water | | Water | | Water | | Water | | Water | |
| Units : | | µg/L | | µg/L | | µg/L | | µg/L | | µg/L | |
| Date Sampled : | | 9/01/2011 | | 9/01/2011 | | 9/01/2011 | | 9/01/2011 | | 9/01/2011 | |
| Time Sampled : | | 11:15 | | 08:40 | | 9:35 | | 10:15 | | 13:20 | |
| Dilution Factor : | | 1.0 | | 1.0 | | 1.0 | | 1.0 | | 1.0 | |
| ANALYTE | RL | Result | Flag | Result | Flag | Result | Flag | Result | Flag | Result | Flag |
| SILVER | 1.0 | | | | | | | 0.058 | J | | |
| ALUMINUM | 30 | 6.8 | B | 14 | B | 3.4 | B | 4200 | | 21 | J |
| ARSENIC | 1.0 | 1.2 | B | 1.7 | B | 1.7 | B | 7.2 | | 1.6 | B |
| BARIUM | 10 | 140 | | 120 | | 260 | | 230 | | 650 | |
| BORON | 5.0 | 46 | J | 27 | J | 35 | J | 30 | J | 42 | J |
| BERYLLIUM | 1.0 | | | | | | | 0.24 | J | | |
| CALCIUM | 100 | 4400 | | 37000 | | 41000 | | 30000 | | 25000 | |
| CADMIUM | 1.0 | | | | | | | | | | |
| COBALT | 0.50 | | | 0.19 | | 0.082 | J | 3.1 | | 0.10 | J |
| CHROMIUM | 2.0 | 2.2 | B | 2.3 | B | 2.3 | B | 3.9 | B | 2.2 | B |
| COPPER | 2.0 | 1.5 | B | 2.5 | B | 7.2 | | 27 | | 1.4 | B |
| IRON | 50 | 250 | | 110 | | 22 | B | 5000 | | 1600 | |
| POTASSIUM | 100 | 1900 | | 1000 | | 1400 | | 2600 | | 1400 | |
| MAGNESIUM | 100 | 1200 | J | 8700 | J | 8700 | J | 7100 | J | 8100 | J |
| MANGANESE | 0.50 | 1.8 | J | 95 | J | 1.4 | J | 200 | J | 72 | J |
| SODIUM | 100 | 16000 | | 10000 | | 10000 | | 8600 | | 12000 | |
| NICKEL | 1.0 | 0.66 | B | 0.50 | B | 0.49 | B | 4.2 | | 0.73 | B |
| LEAD | 1.0 | 0.078 | J | 1.5 | | 0.31 | J | 10 | | 0.25 | J |
| MOLYBDENUM | 5.0 | 0.60 | B | 6.3 | | 1.5 | B | 1.2 | B | 1.6 | B |
| ANTIMONY | 2.0 | 0.094 | B | 0.12 | B | 0.091 | B | 1.9 | B | 0.31 | B |
| SELENIUM | 5.0 | 2.0 | B | 1.5 | B | 1.6 | B | 1.2 | B | 0.71 | B |
| THALLIUM | 1.0 | | | 0.27 | J | 0.12 | J | 0.078 | J | | |
| VANADIUM | 1.0 | | | | | 0.86 | B | 5.0 | | 0.40 | B |
| ZINC | 5.0 | 10 | | 24 | | 7.7 | | 140 | | 8.5 | |
| MERCURY | 0.20 | | | | | | | | | | |

RL = Reporting Limit

SEE NARRATIVE FOR CODE DEFINITIONS

To calculate sample quantitation limits: (RL * Dilution Factor)

Revised 09/99

DATA SUMMARY FORM: INORGANIC

Page 2 of 6

Project #: 180-3644-1

SDG: H-1

Site : DIMOCK

Lab. : TALP

TOTAL METALS

| Sample Number: | | R-2 | | RD-1 | | RU-1 | | KDE-1 | | N-1 | |
|-------------------|------|------------|------|------------|------|------------|------|------------|------|-------------|------|
| Station Location: | | R-2 | | RD-1 | | RU-1 | | KDE-1 | | N-1 | |
| Laboratory ID: | | 180-3644-6 | | 180-3644-7 | | 180-3644-8 | | 180-3644-9 | | 180-3644-10 | |
| Field QC: | | | | | | | | | | | |
| Matrix : | | Water | | Water | | Water | | Water | | Water | |
| Units : | | µg/L | | µg/L | | µg/L | | µg/L | | µg/L | |
| Date Sampled : | | 9/01/2011 | | 9/01/2011 | | 9/01/2011 | | 9/01/2011 | | 9/01/2011 | |
| Time Sampled : | | 13:40 | | 14:10 | | 14:45 | | 15:35 | | 17:45 | |
| Dilution Factor : | | 1.0 | | 1.0 | | 1.0 | | 1.0 | | 1.0 | |
| ANALYTE | RL | Result | Flag | Result | Flag | Result | Flag | Result | Flag | Result | Flag |
| SILVER | 1.0 | | | | | | | | | | |
| ALUMINUM | 30 | 23 | J | 33 | | 32 | | 7.6 | B | 24 | J |
| ARSENIC | 1.0 | 1.1 | B | 0.78 | B | 0.52 | B | 0.67 | B | | |
| BARIUM | 10 | 640 | | 33 | | 33 | | 100 | | 0.12 | J |
| BORON | 5.0 | 41 | J | 14 | J | 15 | J | 180 | J | 8.4 | J |
| BERYLLIUM | 1.0 | | | | | | | | | | |
| CALCIUM | 100 | 26000 | | 9300 | | 9400 | | 1500 | | 11000 | |
| CADMIUM | 1.0 | | | | | | | | | | |
| COBALT | 0.50 | 0.098 | J | | | | | | | | |
| CHROMIUM | 2.0 | 1.2 | B | 1.2 | B | 1.1 | B | 1.3 | B | 1.3 | B |
| COPPER | 2.0 | 2.2 | B | 0.90 | B | 1.0 | B | 4.2 | | 0.33 | B |
| IRON | 50 | 1400 | | 40 | B | 36 | B | 15 | B | 11 | B |
| POTASSIUM | 100 | 1400 | | 1400 | | 1400 | | 1500 | | 70 | B |
| MAGNESIUM | 100 | 8100 | J | 1900 | J | 1900 | J | 160 | J | 3800 | J |
| MANGANESE | 0.50 | 74 | J | 3.1 | J | 20 | J | 2.4 | J | 0.056 | B |
| SODIUM | 100 | 12000 | | 5500 | | 5500 | | 100000 | | 9600 | |
| NICKEL | 1.0 | 0.38 | B | 0.21 | B | 0.25 | B | | | | |
| LEAD | 1.0 | 0.58 | J | 0.055 | J | 0.064 | J | 0.080 | J | | |
| MOLYBDENUM | 5.0 | 1.2 | B | 0.22 | B | | | 0.33 | B | | |
| ANTIMONY | 2.0 | 0.16 | B | 0.14 | B | 0.084 | B | 0.058 | B | 0.19 | B |
| SELENIUM | 5.0 | 1.0 | B | 1.3 | B | 1.2 | B | 0.93 | B | 1.1 | B |
| THALLIUM | 1.0 | | | | | | | | | | |
| VANADIUM | 1.0 | | | 0.44 | B | | | | | | |
| ZINC | 5.0 | 56 | | 1.3 | J | 1.9 | J | 2.2 | J | | |
| MERCURY | 0.20 | | | | | | | | | | |

RL = Reporting Limit

SEE NARRATIVE FOR CODE DEFINITIONS

To calculate sample quantitation limits: (RL * Dilution Factor)

Revised 09/99

DATA SUMMARY FORM: INORGANIC

Page 3 of 6

Project #: 180-3644-1

SDG: H-1

Site : DIMOCK

Lab. : TALP

TOTAL METALS

| Sample Number: | D-1 | | | | | | | | | | |
|-------------------|-------------|--------|------|--------|------|--------|------|--------|------|--------|------|
| Station Location: | D-1 | | | | | | | | | | |
| Laboratory ID: | 180-3644-11 | | | | | | | | | | |
| Field QC: | | | | | | | | | | | |
| Matrix : | Water | | | | | | | | | | |
| Units : | µg/L | | | | | | | | | | |
| Date Sampled : | 9/01/2011 | | | | | | | | | | |
| Time Sampled : | 18:15 | | | | | | | | | | |
| Dilution Factor : | 1.0 | | | | | | | | | | |
| ANALYTE | RL | Result | Flag | Result | Flag | Result | Flag | Result | Flag | Result | Flag |
| SILVER | 1.0 | | | | | | | | | | |
| ALUMINUM | 30 | 2.7 | B | | | | | | | | |
| ARSENIC | 1.0 | 0.73 | B | | | | | | | | |
| BARIUM | 10 | 0.15 | J | | | | | | | | |
| BORON | 5.0 | 5.0 | B | | | | | | | | |
| BERYLLIUM | 1.0 | | | | | | | | | | |
| CALCIUM | 100 | 250 | | | | | | | | | |
| CADMIUM | 1.0 | | | | | | | | | | |
| COBALT | 0.50 | | | | | | | | | | |
| CHROMIUM | 2.0 | 1.3 | B | | | | | | | | |
| COPPER | 2.0 | 0.31 | B | | | | | | | | |
| IRON | 50 | 13 | B | | | | | | | | |
| POTASSIUM | 100 | 3100 | | | | | | | | | |
| MAGNESIUM | 100 | 3000 | J | | | | | | | | |
| MANGANESE | 0.50 | 0.054 | B | | | | | | | | |
| SODIUM | 100 | 1600 | | | | | | | | | |
| NICKEL | 1.0 | | | | | | | | | | |
| LEAD | 1.0 | | | | | | | | | | |
| MOLYBDENUM | 5.0 | | | | | | | | | | |
| ANTIMONY | 2.0 | 0.15 | B | | | | | | | | |
| SELENIUM | 5.0 | 0.60 | B | | | | | | | | |
| THALLIUM | 1.0 | | | | | | | | | | |
| VANADIUM | 1.0 | 0.24 | B | | | | | | | | |
| ZINC | 5.0 | 3.4 | J | | | | | | | | |
| MERCURY | 0.20 | | | | | | | | | | |

RL = Reporting Limit

SEE NARRATIVE FOR CODE DEFINITIONS

To calculate sample quantitation limits: (RL * Dilution Factor)

Revised 09/99

DATA SUMMARY FORM: INORGANIC

Page 4 of 6

Project #: 180-3644-1

SDG: H-1

Site : DIMOCK

Lab. : TALP

DISSOLVED METALS

| Sample Number: | | H-1 | | FH-1 | | FPT-1 | | S-1 | | R-1 | |
|-------------------|------|------------|------|------------|------|------------|------|------------|------|------------|------|
| Station Location: | | H-1 | | FH-1 | | FPT-1 | | S-1 | | R-1 | |
| Laboratory ID: | | 180-3644-1 | | 180-3644-2 | | 180-3644-3 | | 180-3644-4 | | 180-3644-5 | |
| Field QC: | | | | | | | | | | | |
| Matrix : | | Water | | Water | | Water | | Water | | Water | |
| Units : | | µg/L | | µg/L | | µg/L | | µg/L | | µg/L | |
| Date Sampled : | | 9/01/2011 | | 9/01/2011 | | 9/01/2011 | | 9/01/2011 | | 9/01/2011 | |
| Time Sampled : | | 11:15 | | 08:40 | | 9:35 | | 10:15 | | 13:20 | |
| Dilution Factor : | | 1.0 | | 1.0 | | 1.0 | | 1.0 | | 1.0 | |
| ANALYTE | RL | Result | Flag | Result | Flag | Result | Flag | Result | Flag | Result | Flag |
| SILVER | 1.0 | 0.065 | B | 0.042 | B | | | | | | |
| ALUMINUM | 30 | 3.6 | B | 4.1 | B | | | 3.8 | B | | |
| ARSENIC | 1.0 | 0.82 | J | 0.85 | J | | | 1.1 | | 1.0 | |
| BARIUM | 10 | 120 | | 110 | | 230 | | 150 | | 580 | |
| BORON | 5.0 | 36 | | 13 | B | 230 | | 29 | | 38 | |
| BERYLLIUM | 1.0 | | | | | | | 0.059 | B | | |
| CALCIUM | 100 | 3900 | | 35000 | | 37000 | | 29000 | | 23000 | |
| CADMIUM | 1.0 | | | | | | | | | | |
| COBALT | 0.50 | | | 0.039 | B | 0.090 | B | 0.056 | B | 0.097 | B |
| CHROMIUM | 2.0 | 1.9 | B | 1.6 | B | 1.9 | B | 2.0 | B | 1.8 | B |
| COPPER | 2.0 | 1.2 | B | 1.3 | B | 5.6 | | 2.2 | B | 1.1 | B |
| IRON | 50 | 21 | B | 7.4 | B | 50 | B | | | | |
| POTASSIUM | 100 | 1600 | | 980 | | 1500 | | 1100 | | 1300 | |
| MAGNESIUM | 100 | 1000 | | 7700 | | 7400 | | 5900 | | 7000 | |
| MANGANESE | 0.50 | 0.68 | B | 3.4 | | 0.37 | B | 0.19 | B | 68 | |
| SODIUM | 100 | 14000 | | 9600 | | 10000 | | 8700 | | 11000 | |
| NICKEL | 1.0 | | | | | | | 0.19 | J | 0.17 | J |
| LEAD | 1.0 | 0.045 | B | 0.077 | B | 0.098 | B | 0.029 | B | 0.030 | B |
| MOLYBDENUM | 5.0 | 3.7 | B | 1.4 | B | 1.5 | B | 0.90 | B | 1.0 | B |
| ANTIMONY | 2.0 | 0.55 | B | 0.27 | B | 0.30 | B | 0.18 | B | 0.070 | B |
| SELENIUM | 5.0 | 0.53 | J | | | | | | | | |
| THALLIUM | 1.0 | 0.14 | B | 0.073 | B | 0.11 | B | 0.10 | B | 0.043 | B |
| VANADIUM | 1.0 | | | | | 0.21 | B | 0.33 | B | | |
| ZINC | 5.0 | 3.9 | B | 8.3 | B | 13 | | 48 | | 3.4 | B |
| MERCURY | 0.20 | | | | | | | | | | |

RL = Reporting Limit

SEE NARRATIVE FOR CODE DEFINITIONS

To calculate sample quantitation limits: (RL * Dilution Factor)

Revised 09/99

DATA SUMMARY FORM: INORGANIC

Page 5 of 6

Project #: 180-3644-1

SDG: H-1

Site : DIMOCK

Lab. : TALP

DISSOLVED METALS

| Sample Number: | | R-2 | | RD-1 | | RU-1 | | KDE-1 | | N-1 | |
|-------------------|------|------------|------|------------|------|------------|------|------------|------|-------------|------|
| Station Location: | | R-2 | | RD-1 | | RU-1 | | KDE-1 | | N-1 | |
| Laboratory ID: | | 180-3644-6 | | 180-3644-7 | | 180-3644-8 | | 180-3644-9 | | 180-3644-10 | |
| Field QC: | | | | | | | | | | | |
| Matrix : | | Water | | Water | | Water | | Water | | Water | |
| Units : | | µg/L | | µg/L | | µg/L | | µg/L | | µg/L | |
| Date Sampled : | | 9/01/2011 | | 9/01/2011 | | 9/01/2011 | | 9/01/2011 | | 9/01/2011 | |
| Time Sampled : | | 13:40 | | 14:10 | | 14:45 | | 15:35 | | 17:45 | |
| Dilution Factor : | | 1.0 | | 1.0 | | 1.0 | | 1.0 | | 1.0 | |
| ANALYTE | RL | Result | Flag | Result | Flag | Result | Flag | Result | Flag | Result | Flag |
| SILVER | 1.0 | | | | | | | | | | |
| ALUMINUM | 30 | 11 | B | 6.1 | B | 7.4 | B | 4.3 | B | | |
| ARSENIC | 1.0 | 1.2 | | 0.40 | J | 0.78 | J | 0.49 | J | | |
| BARIUM | 10 | 560 | | 30 | | 32 | | 86 | | 0.16 | B |
| BORON | 5.0 | 36 | | 11 | B | 13 | B | 160 | | 7.3 | B |
| BERYLLIUM | 1.0 | | | | | | | | | | |
| CALCIUM | 100 | 23000 | | 8600 | | 8800 | | 1400 | | 11000 | |
| CADMIUM | 1.0 | | | | | | | | | | |
| COBALT | 0.50 | 0.071 | B | 0.032 | B | | | | | | |
| CHROMIUM | 2.0 | 1.6 | B | 1.0 | B | 1.2 | B | 1.2 | B | 1.6 | B |
| COPPER | 2.0 | 0.91 | B | 1.6 | B | 1.7 | B | 2.3 | B | 1.2 | B |
| IRON | 50 | 14 | B | | | | | | | | |
| POTASSIUM | 100 | 1300 | | 1300 | | 1300 | | 1300 | | 200 | |
| MAGNESIUM | 100 | 6800 | | 1700 | | 1700 | | 170 | | 3500 | |
| MANGANESE | 0.50 | 67 | | 0.27 | B | 1.5 | B | 3.0 | | 0.13 | B |
| SODIUM | 100 | 11000 | | 5500 | | 5600 | | 90000 | | 9600 | |
| NICKEL | 1.0 | 0.25 | J | | | | | | | 0.69 | J |
| LEAD | 1.0 | 0.10 | B | 0.029 | B | 0.050 | B | 0.074 | B | 0.036 | B |
| MOLYBDENUM | 5.0 | 0.95 | B | | | | | | | | |
| ANTIMONY | 2.0 | 0.12 | B | 1.2 | B | 0.40 | B | 0.28 | B | 0.30 | B |
| SELENIUM | 5.0 | 0.44 | J | | | 0.43 | J | | | | |
| THALLIUM | 1.0 | 0.050 | B | 0.025 | B | 0.020 | B | | | 0.015 | B |
| VANADIUM | 1.0 | | | | | 0.50 | B | | | | |
| ZINC | 5.0 | 9.1 | | 3.0 | B | 4.8 | B | 3.8 | B | 3.4 | B |
| MERCURY | 0.20 | | | | | | | | | | |

RL = Reporting Limit

SEE NARRATIVE FOR CODE DEFINITIONS

To calculate sample quantitation limits: (RL * Dilution Factor)

Revised 09/99

DATA SUMMARY FORM: INORGANIC

Page 6 of 6

Project #: 180-3644-1

SDG: H-1

Site : DIMOCK

Lab. : TALP

DISSOLVED METALS

| Sample Number: | | D-1 | | | | | | | | | |
|-------------------|------|-------------|------|--------|------|--------|------|--------|------|--------|------|
| Station Location: | | D-1 | | | | | | | | | |
| | | 180-3644-11 | | | | | | | | | |
| Matrix : | | Water | | | | | | | | | |
| Units : | | µg/L | | | | | | | | | |
| Date Sampled : | | 9/01/2011 | | | | | | | | | |
| Time Sampled : | | 18:15 | | | | | | | | | |
| Dilution Factor : | | 1.0 | | | | | | | | | |
| ANALYTE | RL | Result | Flag | Result | Flag | Result | Flag | Result | Flag | Result | Flag |
| SILVER | 1.0 | | | | | | | | | | |
| ALUMINUM | 30 | | | | | | | | | | |
| ARSENIC | 1.0 | | | | | | | | | | |
| BARIUM | 10 | 0.18 | B | | | | | | | | |
| BORON | 5.0 | 3.0 | B | | | | | | | | |
| BERYLLIUM | 1.0 | | | | | | | | | | |
| CALCIUM | 100 | 400 | | | | | | | | | |
| CADMIUM | 1.0 | | | | | | | | | | |
| COBALT | 0.50 | | | | | | | | | | |
| CHROMIUM | 2.0 | 1.5 | B | | | | | | | | |
| COPPER | 2.0 | 1.1 | B | | | | | | | | |
| IRON | 50 | | | | | | | | | | |
| POTASSIUM | 100 | 2900 | | | | | | | | | |
| MAGNESIUM | 100 | 2700 | | | | | | | | | |
| MANGANESE | 0.50 | 0.16 | B | | | | | | | | |
| SODIUM | 100 | 1600 | | | | | | | | | |
| NICKEL | 1.0 | | | | | | | | | | |
| LEAD | 1.0 | 0.037 | B | | | | | | | | |
| MOLYBDENUM | 5.0 | | | | | | | | | | |
| ANTIMONY | 2.0 | 0.21 | B | | | | | | | | |
| SELENIUM | 5.0 | 0.59 | J | | | | | | | | |
| THALLIUM | 1.0 | | | | | | | | | | |
| VANADIUM | 1.0 | | | | | | | | | | |
| ZINC | 5.0 | 3.1 | B | | | | | | | | |
| MERCURY | 0.20 | | | | | | | | | | |

RL = Reporting Limit

SEE NARRATIVE FOR CODE DEFINITIONS

To calculate sample quantitation limits: (RL * Dilution Factor)

Revised 09/99

APPENDIX C

Chain of Custody (COC) Records

Chain of Custody Record

Water

TestAmerica

THE LEADER IN ENVIRONMENTAL TESTING

TestAmerica Laboratory location:

Pittsburgh

Regulatory program:

☐ DW ☐ NPDES ☐ RCRA ☐ Other

TestAmerica Laboratories, Inc.

| | | | | | | | | | |
|-----------------------------------|--|-----------------------------|--|--|--|---|--|--|--|
| Client Contact | | Client Project Manager: | | Site Contact: | | Lab Contact: | | COC No: | |
| URS Corporation | | Thomas Merski | | Jim Pinta | | Carrie Gamber | | 005410 | |
| Address: | | Telephone: | | Telephone: | | Telephone: | | | |
| 501 Holiday Dr., Suite 300 | | 412-503-4603 | | 412-860-6342 (cell) | | 412-963-7058 | | | |
| City/State/Zip: | | Email: | | Email: | | Email: | | | |
| Pittsburgh, PA 15220 | | James.Pinta@urscorp.com | | David.Blanchild@urscorp.com | | | | | |
| Phone: | | Method of Shipment/Carrier: | | Shipping/Tracking No: | | | | | |
| 412-503-4700 | | FedEx Next Day Air | | | | | | | |
| Project Name: | | Project Number: | | TAT if different than below | | Analysis: | | | |
| K&L Gates Focused Site Assessment | | 39938688.00001 | | <input type="checkbox"/> 3 weeks <input type="checkbox"/> 2 weeks <input type="checkbox"/> 1 week <input type="checkbox"/> 2 days <input type="checkbox"/> 1 day | | VOCs SVOCs Metals (not dissolved) PCBs Pesticides Herbicides Fungicides PAHs GYCOLS DROGRO MBAS | | | |
| PO# | | Sample Identification | | Sample Date | | Sample Time | | Sample Specific Notes / Special Instructions | |
| 644872 | | | | | | | | | |
| | | H-1 | | 9/1/11 | | 11:55 AM | | X X X X X X X X X X | |
| | | FH-1 | | 9/1/11 | | 8:40 AM | | X X X X X X X X X X | |
| | | FPT-1 | | 9/1/11 | | 7:35 AM | | X X X X X X X X X X | |
| | | S-1 | | 9/1/11 | | 10:15 AM | | X X X X X X X X X X | |
| | | R-1 | | 9/1/11 | | 1:20 PM | | X X X X X X X X X X | |
| | | R-2 | | 9/1/11 | | 1:40 PM | | X X X X X X X X X X | |
| | | RD-1 | | 9/1/11 | | 2:10 PM | | X X X X X X X X X X | |
| | | RHE-1 | | 9/1/11 | | 2:45 PM | | X X X X X X X X X X | |
| | | N-1 | | 9/1/11 | | 5:45 PM | | X X X X X X X X X X | |
| | | D-1 | | 9/1/11 | | 8:15 PM | | X X X X X X X X X X | |
| Possible Hazard Identification | | Non-Hazard | | Flammable | | Skin Irritant | | Poison B | |
| | | <input type="checkbox"/> | | <input type="checkbox"/> | | <input type="checkbox"/> | | <input type="checkbox"/> | |
| Special Instructions/Comments | | Unknown | | Sample Disposal (A fee may be assessed if samples are retained longer than 1 month) | | Return to Client | | Dispose By Lab | |
| | | <input type="checkbox"/> | | <input type="checkbox"/> | | <input type="checkbox"/> | | <input type="checkbox"/> | |
| Relinquished by: | | Company: | | Date/Time: | | Received by: | | Date/Time: | |
| [Signature] | | URS Corp | | 9/1/11 | | [Signature] | | 9/6/11 0900 | |
| Relinquished by: | | Company: | | Date/Time: | | Received by: | | Date/Time: | |
| [Signature] | | URS Corp | | 9/6/11 | | [Signature] | | 9/6/11 0900 | |
| Relinquished by: | | Company: | | Date/Time: | | Received in Laboratory by: | | Date/Time: | |
| [Signature] | | URS Corp | | 9/6/11 | | [Signature] | | 9/6/11 0900 | |

©2008 TestAmerica Laboratories, Inc. All rights reserved. TestAmerica & Chain of Custody are trademarks of TestAmerica Laboratories, Inc.

TAL 0018-1 (04/10)

CABOT-EPA 003754

Page 314 of 317

9/6/2011

Login Container Summary Report

180-3644

Temperature readings: _____

| <u>Client Sample ID</u> | <u>Lab ID</u> | <u>Container Type</u> | <u>Container</u> <u>pH</u> | <u>Preservative</u> <u>Added (mls)</u> | <u>Lot #</u> |
|-------------------------|---------------|-------------------------------------|-------------------------------|---|--------------|
| H-1 | 180-3644-A-1 | Plastic 1 liter - unpreserved | _____ | _____ | _____ |
| H-1 | 180-3644-B-1 | Amber Glass 1 liter - Hydrochloric | 2 | _____ | _____ |
| H-1 | 180-3644-C-1 | Amber Glass 1 liter - Hydrochloric | 2 | _____ | _____ |
| H-1 | 180-3644-D-1 | Amber Glass 1 liter - Sulfuric Acid | 2 | _____ | _____ |
| H-1 | 180-3644-E-1 | Amber Glass 1 liter - Sulfuric Acid | 2 | _____ | _____ |
| H-1 | 180-3644-F-1 | Amber Glass 1 liter - unpreserved | _____ | _____ | _____ |
| H-1 | 180-3644-G-1 | Amber Glass 1 liter - unpreserved | _____ | _____ | _____ |
| H-1 | 180-3644-H-1 | Plastic 500ml - with Nitric Acid | 2 | _____ | _____ |
| H-1 | 180-3644-I-1 | Plastic 500ml - unpreserved | _____ | _____ | _____ |
| H-1 | 180-3644-J-1 | Plastic 250ml - with Sulfuric Acid | 2 | _____ | _____ |
| H-1 | 180-3644-K-1 | Voa Vial 40ml - unpreserved | _____ | _____ | _____ |
| H-1 | 180-3644-L-1 | Voa Vial 40ml - unpreserved | _____ | _____ | _____ |
| H-1 | 180-3644-M-1 | Voa Vial 40ml - unpreserved | _____ | _____ | _____ |
| H-1 | 180-3644-N-1 | Voa Vial 40ml - Hydrochloric Acid | p | _____ | _____ |
| H-1 | 180-3644-O-1 | Voa Vial 40ml - Hydrochloric Acid | ↓ | _____ | _____ |
| H-1 | 180-3644-P-1 | Voa Vial 40ml - Hydrochloric Acid | ↓ | _____ | _____ |
| H-1 | 180-3644-Q-1 | Voa Vial 40ml - Hydrochloric Acid | ↓ | _____ | _____ |
| H-1 | 180-3644-R-1 | Voa Vial 40ml - Hydrochloric Acid | ↓ | _____ | _____ |
| H-1 | 180-3644-S-1 | Voa Vial 40ml - Hydrochloric Acid | ↓ | _____ | _____ |
| H-1 | 180-3644-T-1 | Voa Vial 40ml - with Sodium | p | _____ | _____ |
| H-1 | 180-3644-U-1 | Voa Vial 40ml - with Sodium | ↓ | _____ | _____ |
| H-1 | 180-3644-V-1 | Voa Vial 40ml - with Sulfuric Acid | ↓ | _____ | _____ |
| H-1 | 180-3644-W-1 | Voa Vial 40ml - with Sulfuric Acid | ↓ | _____ | _____ |
| FH-1 | 180-3644-A-2 | Plastic 1 liter - unpreserved | _____ | _____ | _____ |
| FH-1 | 180-3644-B-2 | Amber Glass 1 liter - Hydrochloric | 2 | _____ | _____ |
| FH-1 | 180-3644-C-2 | Amber Glass 1 liter - Hydrochloric | 2 | _____ | _____ |
| FH-1 | 180-3644-D-2 | Amber Glass 1 liter - Sulfuric Acid | 2 | _____ | _____ |
| FH-1 | 180-3644-E-2 | Amber Glass 1 liter - Sulfuric Acid | 2 | _____ | _____ |
| FH-1 | 180-3644-F-2 | Amber Glass 1 liter - unpreserved | _____ | _____ | _____ |
| FH-1 | 180-3644-G-2 | Amber Glass 1 liter - unpreserved | _____ | _____ | _____ |
| FH-1 | 180-3644-H-2 | Plastic 500ml - with Nitric Acid | 2 | _____ | _____ |
| FH-1 | 180-3644-I-2 | Plastic 500ml - unpreserved | _____ | _____ | _____ |
| FH-1 | 180-3644-J-2 | Plastic 250ml - with Sulfuric Acid | 2 | _____ | _____ |
| FH-1 | 180-3644-K-2 | Voa Vial 40ml - unpreserved | _____ | _____ | _____ |
| FH-1 | 180-3644-L-2 | Voa Vial 40ml - unpreserved | _____ | _____ | _____ |
| FH-1 | 180-3644-M-2 | Voa Vial 40ml - unpreserved | _____ | _____ | _____ |

| <u>Client Sample ID</u> | <u>Lab ID</u> | <u>Container Type</u> | <u>Container</u> <u>pH</u> | <u>Preservative</u> <u>Added (mls)</u> | <u>Lot #</u> |
|-------------------------|---------------|---|-------------------------------|---|--------------|
| FH-1 | 180-3644-N-2 | Voa Vial 40ml - Hydrochloric Acid | P | | |
| FH-1 | 180-3644-O-2 | Voa Vial 40ml - Hydrochloric Acid | | | |
| FH-1 | 180-3644-P-2 | Voa Vial 40ml - Hydrochloric Acid | | | |
| FH-1 | 180-3644-Q-2 | Voa Vial 40ml - Hydrochloric Acid | | | |
| FH-1 | 180-3644-R-2 | Voa Vial 40ml - Hydrochloric Acid | | | |
| FH-1 | 180-3644-S-2 | Voa Vial 40ml - Hydrochloric Acid | | | |
| FH-1 | 180-3644-T-2 | Voa Vial 40ml - with Sodium | | | |
| FH-1 | 180-3644-U-2 | Voa Vial 40ml - with Sodium | | | |
| FH-1 | 180-3644-V-2 | Voa Vial 40ml - with Sulfuric Acid | | | |
| FH-1 | 180-3644-W-2 | Voa Vial 40ml - with Sulfuric Acid | ↓ | | |
| FPT-1 | 180-3644-A-3 | Plastic 1 liter - unpreserved | | | |
| FPT-1 | 180-3644-B-3 | Amber Glass 1 liter - Hydrochloric | 2 | | |
| FPT-1 | 180-3644-C-3 | Amber Glass 1 liter - Hydrochloric | 2 | | |
| FPT-1 | 180-3644-D-3 | Amber Glass 1 liter - Sulfuric Acid | 2 | | |
| FPT-1 | 180-3644-F-3 | Amber Glass 1 liter - unpreserved <i>sulfuric acid</i> | 2 | | |
| FPT-1 | 180-3644-G-3 | Amber Glass 1 liter - unpreserved | | | |
| FPT-1 | 180-3644-H-3 | Plastic 500ml - with Nitric Acid | 2 | | |
| FPT-1 | 180-3644-I-3 | Plastic 500ml - unpreserved | | | |
| FPT-1 | 180-3644-J-3 | Plastic 250ml - with Sulfuric Acid | 2 | | |
| FPT-1 | 180-3644-K-3 | Voa Vial 40ml - unpreserved | | | |
| FPT-1 | 180-3644-L-3 | Voa Vial 40ml - unpreserved | | | |
| FPT-1 | 180-3644-M-3 | Voa Vial 40ml - unpreserved | | | |
| FPT-1 | 180-3644-N-3 | Voa Vial 40ml - Hydrochloric Acid | P | | |
| FPT-1 | 180-3644-O-3 | Voa Vial 40ml - Hydrochloric Acid | | | |
| FPT-1 | 180-3644-P-3 | Voa Vial 40ml - Hydrochloric Acid | | | |
| FPT-1 | 180-3644-Q-3 | Voa Vial 40ml - Hydrochloric Acid | | | |
| FPT-1 | 180-3644-R-3 | Voa Vial 40ml - Hydrochloric Acid | | | |
| FPT-1 | 180-3644-S-3 | Voa Vial 40ml - Hydrochloric Acid | | | |
| FPT-1 | 180-3644-T-3 | Voa Vial 40ml - with Sodium | | | |
| FPT-1 | 180-3644-U-3 | Voa Vial 40ml - with Sodium | | | |
| FPT-1 | 180-3644-V-3 | Voa Vial 40ml - with Sulfuric Acid | | | |
| FPT-1 | 180-3644-W-3 | Voa Vial 40ml - with Sulfuric Acid | ✓ | | |
| S-1 | 180-3644-A-4 | Plastic 1 liter - unpreserved | | | |
| S-1 | 180-3644-B-4 | Amber Glass 1 liter - Hydrochloric | 2 | | |
| S-1 | 180-3644-C-4 | Amber Glass 1 liter - Hydrochloric | 2 | | |
| S-1 | 180-3644-D-4 | Amber Glass 1 liter - Sulfuric Acid | 2 | | |
| S-1 | 180-3644-E-4 | Amber Glass 1 liter - Sulfuric Acid | 2 | | |
| S-1 | 180-3644-F-4 | Amber Glass 1 liter - unpreserved | | | |
| S-1 | 180-3644-G-4 | Amber Glass 1 liter - unpreserved | | | |

| <u>Client Sample ID</u> | <u>Lab ID</u> | <u>Container Type</u> | <u>Container</u> <u>pH</u> | <u>Preservative</u> <u>Added (mls)</u> | <u>Lot #</u> |
|-------------------------|---------------|------------------------------------|-------------------------------|---|--------------|
| S-1 | 180-3644-H-4 | Plastic 500ml - with Nitric Acid | ~ | | |
| S-1 | 180-3644-I-4 | Plastic 500ml - unpreserved | | | |
| S-1 | 180-3644-J-4 | Plastic 250ml - with Sulfuric Acid | ~ | | |
| S-1 | 180-3644-K-4 | Voa Vial 40ml - unpreserved | | | |
| S-1 | 180-3644-L-4 | Voa Vial 40ml - unpreserved | | | |
| S-1 | 180-3644-M-4 | Voa Vial 40ml - unpreserved | | | |
| S-1 | 180-3644-N-4 | Voa Vial 40ml - Hydrochloric Acid | p | | |
| S-1 | 180-3644-O-4 | Voa Vial 40ml - Hydrochloric Acid | | | |
| S-1 | 180-3644-P-4 | Voa Vial 40ml - Hydrochloric Acid | | | |
| S-1 | 180-3644-Q-4 | Voa Vial 40ml - Hydrochloric Acid | | | |
| S-1 | 180-3644-R-4 | Voa Vial 40ml - Hydrochloric Acid | | | |
| S-1 | 180-3644-S-4 | Voa Vial 40ml - Hydrochloric Acid | | | |
| S-1 | 180-3644-T-4 | Voa Vial 40ml - with Sodium | | | |
| S-1 | 180-3644-U-4 | Voa Vial 40ml - with Sodium | | | |
| S-1 | 180-3644-V-4 | Voa Vial 40ml - with Sulfuric Acid | | | |
| S-1 | 180-3644-W-4 | Voa Vial 40ml - with Sulfuric Acid | ↓ | | |
| R-1 | 180-3644-A-5 | Plastic 1 liter - unpreserved | | | |
| R-1 | 180-3644-B-5 | Amber Glass 1 liter - Hydrochloric | ~ | | |
| R-1 | 180-3644-C-5 | Amber Glass 1 liter - Hydrochloric | ~ | | |
| R-1 | 180-3644-D-5 | Amber Glass 1 liter - unpreserved | | | |
| R-1 | 180-3644-E-5 | Amber Glass 1 liter - unpreserved | | | |
| R-1 | 180-3644-F-5 | Amber Glass 1 liter - unpreserved | | | |
| R-1 | 180-3644-G-5 | Amber Glass 1 liter - unpreserved | | | |
| R-1 | 180-3644-H-5 | Plastic 500ml - with Nitric Acid | ~ | | |
| R-1 | 180-3644-I-5 | Plastic 500ml - unpreserved | | | |
| R-1 | 180-3644-J-5 | Plastic 250ml - with Sulfuric Acid | ~ | | |
| R-1 | 180-3644-K-5 | Voa Vial 40ml - unpreserved | | | |
| R-1 | 180-3644-L-5 | Voa Vial 40ml - unpreserved | | | |
| R-1 | 180-3644-M-5 | Voa Vial 40ml - unpreserved | | | |
| R-1 | 180-3644-N-5 | Voa Vial 40ml - Hydrochloric Acid | p | | |
| R-1 | 180-3644-O-5 | Voa Vial 40ml - Hydrochloric Acid | | | |
| R-1 | 180-3644-P-5 | Voa Vial 40ml - Hydrochloric Acid | | | |
| R-1 | 180-3644-Q-5 | Voa Vial 40ml - Hydrochloric Acid | | | |
| R-1 | 180-3644-R-5 | Voa Vial 40ml - Hydrochloric Acid | | | |
| R-1 | 180-3644-S-5 | Voa Vial 40ml - Hydrochloric Acid | | | |
| R-1 | 180-3644-T-5 | Voa Vial 40ml - Hydrochloric Acid | | | |
| R-1 | 180-3644-U-5 | Voa Vial 40ml - Hydrochloric Acid | | | |
| R-1 | 180-3644-V-5 | Voa Vial 40ml - with Sulfuric Acid | | | |
| R-1 | 180-3644-W-5 | Voa Vial 40ml - with Sulfuric Acid | ↓ | | |

| <u>Client Sample ID</u> | <u>Lab ID</u> | <u>Container Type</u> | <u>Container</u> <u>pH</u> | <u>Preservative</u> <u>Added (mls)</u> | <u>Lot #</u> |
|-------------------------|---------------|-------------------------------------|-------------------------------|---|--------------|
| R-2 | 180-3644-A-6 | Plastic 1 liter - unpreserved | — | — | — |
| R-2 | 180-3644-B-6 | Amber Glass 1 liter - Hydrochloric | 2 | — | — |
| R-2 | 180-3644-C-6 | Amber Glass 1 liter - Hydrochloric | 2 | — | — |
| R-2 | 180-3644-D-6 | Amber Glass 1 liter - Sulfuric Acid | 2 | — | — |
| R-2 | 180-3644-E-6 | Amber Glass 1 liter - Sulfuric Acid | 2 | — | — |
| R-2 | 180-3644-F-6 | Amber Glass 1 liter - unpreserved | — | — | — |
| R-2 | 180-3644-G-6 | Amber Glass 1 liter - unpreserved | — | — | — |
| R-2 | 180-3644-H-6 | Plastic 500ml - with Nitric Acid | 2 | — | — |
| R-2 | 180-3644-I-6 | Plastic 500ml - unpreserved | — | — | — |
| R-2 | 180-3644-J-6 | Plastic 250ml - with Sulfuric Acid | 2 | — | — |
| R-2 | 180-3644-K-6 | Voa Vial 40ml - unpreserved | — | — | — |
| R-2 | 180-3644-L-6 | Voa Vial 40ml - unpreserved | — | — | — |
| R-2 | 180-3644-M-6 | Voa Vial 40ml - unpreserved | — | — | — |
| R-2 | 180-3644-N-6 | Voa Vial 40ml - Hydrochloric Acid | P | — | — |
| R-2 | 180-3644-O-6 | Voa Vial 40ml - Hydrochloric Acid | | — | — |
| R-2 | 180-3644-P-6 | Voa Vial 40ml - Hydrochloric Acid | | — | — |
| R-2 | 180-3644-Q-6 | Voa Vial 40ml - Hydrochloric Acid | | — | — |
| R-2 | 180-3644-R-6 | Voa Vial 40ml - Hydrochloric Acid | | — | — |
| R-2 | 180-3644-S-6 | Voa Vial 40ml - Hydrochloric Acid | | — | — |
| R-2 | 180-3644-T-6 | Voa Vial 40ml - Hydrochloric Acid | | — | — |
| R-2 | 180-3644-U-6 | Voa Vial 40ml - Hydrochloric Acid | | — | — |
| R-2 | 180-3644-V-6 | Voa Vial 40ml - with Sulfuric Acid | | — | — |
| R-2 | 180-3644-W-6 | Voa Vial 40ml - with Sulfuric Acid | ↓ | — | — |
| RD-1 | 180-3644-A-7 | Plastic 1 liter - unpreserved | — | — | — |
| RD-1 | 180-3644-B-7 | Amber Glass 1 liter - Hydrochloric | 2 | — | — |
| RD-1 | 180-3644-D-7 | Amber Glass 1 liter - Sulfuric Acid | 2 | — | — |
| RD-1 | 180-3644-E-7 | Amber Glass 1 liter - Sulfuric Acid | 2 | — | — |
| RD-1 | 180-3644-F-7 | Amber Glass 1 liter - unpreserved | — | — | — |
| RD-1 | 180-3644-G-7 | Amber Glass 1 liter - unpreserved | — | — | — |
| RD-1 | 180-3644-H-7 | Plastic 500ml - with Nitric Acid | 2 | — | — |
| RD-1 | 180-3644-I-7 | Plastic 500ml - unpreserved | — | — | — |
| RD-1 | 180-3644-J-7 | Plastic 250ml - with Sulfuric Acid | 2 | — | — |
| RD-1 | 180-3644-K-7 | Voa Vial 40ml - unpreserved | — | — | — |
| RD-1 | 180-3644-L-7 | Voa Vial 40ml - unpreserved | — | — | — |
| RD-1 | 180-3644-M-7 | Voa Vial 40ml - unpreserved | — | — | — |
| RD-1 | 180-3644-N-7 | Voa Vial 40ml - Hydrochloric Acid | P | — | — |
| RD-1 | 180-3644-O-7 | Voa Vial 40ml - Hydrochloric Acid | | — | — |
| RD-1 | 180-3644-P-7 | Voa Vial 40ml - Hydrochloric Acid | | — | — |
| RD-1 | 180-3644-Q-7 | Voa Vial 40ml - Hydrochloric Acid | ↓ | — | — |

| <u>Client Sample ID</u> | <u>Lab ID</u> | <u>Container Type</u> | <u>Container</u> <u>pH</u> | <u>Preservative</u> <u>Added (mls)</u> | <u>Lot #</u> |
|-------------------------|---------------|-------------------------------------|-------------------------------|---|--------------|
| RD-1 | 180-3644-R-7 | Voa Vial 40ml - Hydrochloric Acid | P | | |
| RD-1 | 180-3644-S-7 | Voa Vial 40ml - Hydrochloric Acid | | | |
| RD-1 | 180-3644-T-7 | Voa Vial 40ml - with Sodium | | | |
| RD-1 | 180-3644-U-7 | Voa Vial 40ml - with Sodium | | | |
| RD-1 | 180-3644-V-7 | Voa Vial 40ml - with Sulfuric Acid | | | |
| RD-1 | 180-3644-W-7 | Voa Vial 40ml - with Sulfuric Acid | ↓ | | |
| RV-1 | 180-3644-A-8 | Plastic 1 liter - unpreserved | | | |
| RV-1 | 180-3644-B-8 | Amber Glass 1 liter - Hydrochloric | 2 | | |
| RV-1 | 180-3644-C-8 | Amber Glass 1 liter - Hydrochloric | 2 | | |
| RV-1 | 180-3644-D-8 | Amber Glass 1 liter - Sulfuric Acid | 2 | | |
| RV-1 | 180-3644-E-8 | Amber Glass 1 liter - Sulfuric Acid | 2 | | |
| RV-1 | 180-3644-F-8 | Amber Glass 1 liter - unpreserved | 26 | | |
| RV-1 | 180-3644-G-8 | Amber Glass 1 liter - unpreserved | | | |
| RV-1 | 180-3644-H-8 | Plastic 500ml - with Nitric Acid | 2 | | |
| RV-1 | 180-3644-I-8 | Plastic 500ml - unpreserved | | | |
| RV-1 | 180-3644-J-8 | Plastic 250ml - with Sulfuric Acid | 2 | | |
| RV-1 | 180-3644-K-8 | Voa Vial 40ml - unpreserved | | | |
| RV-1 | 180-3644-L-8 | Voa Vial 40ml - unpreserved | | | |
| RV-1 | 180-3644-M-8 | Voa Vial 40ml - unpreserved | | | |
| RV-1 | 180-3644-N-8 | Voa Vial 40ml - Hydrochloric Acid | P | | |
| RV-1 | 180-3644-O-8 | Voa Vial 40ml - Hydrochloric Acid | | | |
| RV-1 | 180-3644-P-8 | Voa Vial 40ml - Hydrochloric Acid | | | |
| RV-1 | 180-3644-Q-8 | Voa Vial 40ml - Hydrochloric Acid | | | |
| RV-1 | 180-3644-R-8 | Voa Vial 40ml - Hydrochloric Acid | | | |
| RV-1 | 180-3644-S-8 | Voa Vial 40ml - Hydrochloric Acid | | | |
| RV-1 | 180-3644-T-8 | Voa Vial 40ml - with Sodium | | | |
| RV-1 | 180-3644-U-8 | Voa Vial 40ml - with Sodium | | | |
| RV-1 | 180-3644-V-8 | Voa Vial 40ml - with Sulfuric Acid | | | |
| RV-1 | 180-3644-W-8 | Voa Vial 40ml - with Sulfuric Acid | ↓ | | |
| KDE-1 | 180-3644-A-9 | Plastic 1 liter - unpreserved | | | |
| KDE-1 | 180-3644-B-9 | Amber Glass 1 liter - Hydrochloric | 2 | | |
| KDE-1 | 180-3644-C-9 | Amber Glass 1 liter - Hydrochloric | 2 | | |
| KDE-1 | 180-3644-D-9 | Amber Glass 1 liter - unpreserved | | | |
| KDE-1 | 180-3644-E-9 | Amber Glass 1 liter - unpreserved | | | |
| KDE-1 | 180-3644-F-9 | Amber Glass 1 liter - unpreserved | | | |
| KDE-1 | 180-3644-G-9 | Amber Glass 1 liter - unpreserved | | | |
| KDE-1 | 180-3644-H-9 | Plastic 500ml - with Nitric Acid | 2 | | |
| KDE-1 | 180-3644-I-9 | Plastic 500ml - unpreserved | | | |
| KDE-1 | 180-3644-J-9 | Plastic 250ml - with Sulfuric Acid | 2 | | |

| <u>Client Sample ID</u> | <u>Lab ID</u> | <u>Container Type</u> | <u>Container</u> pH | <u>Preservative</u> Added (mls) | <u>Lot #</u> |
|-------------------------|---------------|------------------------------------|------------------------|------------------------------------|--------------|
| KDE-1 | 180-3644-K-9 | Voa Vial 40ml - unpreserved | — | — | — |
| KDE-1 | 180-3644-L-9 | Voa Vial 40ml - unpreserved | — | — | — |
| KDE-1 | 180-3644-M-9 | Voa Vial 40ml - unpreserved | — | — | — |
| KDE-1 | 180-3644-N-9 | Voa Vial 40ml - Hydrochloric Acid | P | — | — |
| KDE-1 | 180-3644-O-9 | Voa Vial 40ml - Hydrochloric Acid | | — | — |
| KDE-1 | 180-3644-P-9 | Voa Vial 40ml - Hydrochloric Acid | | — | — |
| KDE-1 | 180-3644-Q-9 | Voa Vial 40ml - Hydrochloric Acid | | — | — |
| KDE-1 | 180-3644-R-9 | Voa Vial 40ml - Hydrochloric Acid | | — | — |
| KDE-1 | 180-3644-S-9 | Voa Vial 40ml - Hydrochloric Acid | | — | — |
| KDE-1 | 180-3644-T-9 | Voa Vial 40ml - Hydrochloric Acid | | — | — |
| KDE-1 | 180-3644-U-9 | Voa Vial 40ml - Hydrochloric Acid | | — | — |
| KDE-1 | 180-3644-V-9 | Voa Vial 40ml - with Sulfuric Acid | ↓ | — | — |
| KDE-1 | 180-3644-W-9 | Voa Vial 40ml - with Sulfuric Acid | ↓ | — | — |
| N-1 | 180-3644-A-10 | Plastic 1 liter - unpreserved | — | — | — |
| N-1 | 180-3644-B-10 | Amber Glass 1 liter - Hydrochloric | 2 | — | — |
| N-1 | 180-3644-C-10 | Amber Glass 1 liter - Hydrochloric | 2 | — | — |
| N-1 | 180-3644-D-10 | Amber Glass 1 liter - unpreserved | — | — | — |
| N-1 | 180-3644-E-10 | Amber Glass 1 liter - unpreserved | — | — | — |
| N-1 | 180-3644-F-10 | Amber Glass 1 liter - unpreserved | — | — | — |
| N-1 | 180-3644-G-10 | Amber Glass 1 liter - unpreserved | — | — | — |
| N-1 | 180-3644-H-10 | Plastic 500ml - with Nitric Acid | 2 | — | — |
| N-1 | 180-3644-I-10 | Plastic 500ml - unpreserved | 2 | — | — |
| N-1 | 180-3644-J-10 | Plastic 250ml - with Sulfuric Acid | — | — | — |
| N-1 | 180-3644-K-10 | Voa Vial 40ml - unpreserved | — | — | — |
| N-1 | 180-3644-L-10 | Voa Vial 40ml - unpreserved | — | — | — |
| N-1 | 180-3644-M-10 | Voa Vial 40ml - unpreserved | — | — | — |
| N-1 | 180-3644-N-10 | Voa Vial 40ml - Hydrochloric Acid | P | — | — |
| N-1 | 180-3644-O-10 | Voa Vial 40ml - Hydrochloric Acid | | — | — |
| N-1 | 180-3644-P-10 | Voa Vial 40ml - Hydrochloric Acid | | — | — |
| N-1 | 180-3644-Q-10 | Voa Vial 40ml - Hydrochloric Acid | | — | — |
| N-1 | 180-3644-R-10 | Voa Vial 40ml - Hydrochloric Acid | | — | — |
| N-1 | 180-3644-S-10 | Voa Vial 40ml - Hydrochloric Acid | | — | — |
| N-1 | 180-3644-T-10 | Voa Vial 40ml - with Sodium | — | 01534401V | — |
| N-1 | 180-3644-U-10 | Voa Vial 40ml - with Sodium | — | 01534401W | — |
| N-1 | 180-3644-V-10 | Voa Vial 40ml - with Sulfuric Acid | ↓ | — | — |
| N-1 | 180-3644-W-10 | Voa Vial 40ml - with Sulfuric Acid | ↓ | — | — |
| D-1 | 180-3644-A-11 | Plastic 1 liter - unpreserved | — | — | — |
| D-1 | 180-3644-B-11 | Amber Glass 1 liter - Hydrochloric | 2 | — | — |
| D-1 | 180-3644-C-11 | Amber Glass 1 liter - Hydrochloric | 2 | — | — |

| <u>Client Sample ID</u> | <u>Lab ID</u> | <u>Container Type</u> | <u>Container</u> <u>pH</u> | <u>Preservative</u> <u>Added (mls)</u> | <u>Lot #</u> |
|-------------------------|---------------|-------------------------------------|-------------------------------|---|--------------|
| D-1 | 180-3644-D-11 | Amber Glass 1 liter - Sulfuric Acid | 2 | | |
| D-1 | 180-3644-E-11 | Amber Glass 1 liter - unpreserved | | | |
| D-1 | 180-3644-F-11 | Amber Glass 1 liter - unpreserved | | | |
| D-1 | 180-3644-G-11 | Amber Glass 1 liter - unpreserved | | | |
| D-1 | 180-3644-H-11 | Plastic 500ml - with Nitric Acid | 2 | | |
| D-1 | 180-3644-I-11 | Plastic 500ml - unpreserved | | | |
| D-1 | 180-3644-J-11 | Plastic 250ml - with Sulfuric Acid | 2 | | |
| D-1 | 180-3644-K-11 | Voa Vial 40ml - unpreserved | | | |
| D-1 | 180-3644-L-11 | Voa Vial 40ml - unpreserved | | | |
| D-1 | 180-3644-M-11 | Voa Vial 40ml - unpreserved | | | |
| D-1 | 180-3644-N-11 | Voa Vial 40ml - Hydrochloric Acid | P | | |
| D-1 | 180-3644-O-11 | Voa Vial 40ml - Hydrochloric Acid | | | |
| D-1 | 180-3644-P-11 | Voa Vial 40ml - Hydrochloric Acid | | | |
| D-1 | 180-3644-Q-11 | Voa Vial 40ml - Hydrochloric Acid | | | |
| D-1 | 180-3644-R-11 | Voa Vial 40ml - Hydrochloric Acid | | | |
| D-1 | 180-3644-S-11 | Voa Vial 40ml - Hydrochloric Acid | | | |
| D-1 | 180-3644-T-11 | Voa Vial 40ml - with Sodium | | 0133440W | |
| D-1 | 180-3644-U-11 | Voa Vial 40ml - with Sodium | | 0133440W | |
| D-1 | 180-3644-V-11 | Voa Vial 40ml - with Sulfuric Acid | | 58032-11 | |
| D-1 | 180-3644-W-11 | Voa Vial 40ml - with Sulfuric Acid | | 58032-11 | |
| TRIP BLANK | 180-3644-A-12 | Voa Vial 40ml - Hydrochloric Acid | | | |
| TRIP BLANK | 180-3644-B-12 | Voa Vial 40ml - Hydrochloric Acid | | | |
| TRIP BLANK | 180-3644-C-12 | Voa Vial 40ml - Hydrochloric Acid | 2 | | |

APPENDIX D

Laboratory Case Narrative



THE LEADER IN ENVIRONMENTAL TESTING

ANALYTICAL REPORT

Job Number: 180-3644-1

Job Description: Focused Site Assessment

For:

URS Corporation

Foster Plaza 4

501 Holiday Drive, Suite 300

Pittsburgh, PA 15220

Attention: Mr. James Pinta, Jr.

Approved for release:
Carrie L. Gamber
Project Manager II
10/5/2011 2:08 PM

Carrie L. Gamber

Project Manager II

carrie.gamber@testamericainc.com

10/05/2011

The test results in this report meet all NELAP requirements for parameters for which accreditation is required or available. Any exceptions to the NELAP requirements are noted in this report. Pursuant to NELAP, this report may not be reproduced, except in full, without the written approval of the laboratory. This report is confidential and is intended for the sole use of TestAmerica and its client. All questions regarding this report should be directed to the TestAmerica Project Manager or designee who has signed this report.

TestAmerica Laboratories, Inc.

TestAmerica Pittsburgh 301 Alpha Drive, RIDC Park, Pittsburgh, PA 15238

Tel (412) 963-7058 Fax (412) 963-2468 www.testamericainc.com

CASE NARRATIVE

Client: URS Corporation

Project: Focused Site Assessment

Report Number: 180-3644-1

With the exceptions noted as flags or footnotes, standard analytical protocols were followed in the analysis of the samples and no problems were encountered or anomalies observed. In addition all laboratory quality control samples were within established control limits, with any exceptions noted below. Each sample was analyzed to achieve the lowest possible reporting limit within the constraints of the method. In some cases, due to interference or analytes present at high concentrations, samples were diluted. For diluted samples, the reporting limits are adjusted relative to the dilution required.

Calculations are performed before rounding to avoid round-off errors in calculated results.

All holding times were met and proper preservation noted for the methods performed on these samples, unless otherwise detailed in the individual sections below.

RECEIPT

The samples were received on 09/06/2011; the samples arrived in good condition, properly preserved and on ice.

The laboratory received a broken 1-liter bottle for sample FPT-1 (180-3644-3).

The laboratory did not receive sulfuric acid preserved vials for samples FPT-1 (180-3644-3) and D-1 (180-3644-11). An aliquot was taken from an un-preserved 1-liter amber glass bottle and preserved for TOC.

No Sodium thiosulfate vials were received for samples FPT-1 (180-3644-3), D-1 (180-3644-11), and N-1 (180-3644-10). An aliquot was taken from an un-preserved 1-liter amber glass bottle and preserved for method 8011.

The DRO volume for sample FPT-1 (180-3644-3) was received in a 1-liter plastic container.

LOW LEVEL VOLATILE ORGANIC COMPOUNDS

Chloroform, Naphthalene and Toluene were detected in method blank MB 180-14017/3 at levels that were above the method detection limit but below the reporting limit. The values should be considered estimates, and have been flagged "J". If the associated sample reported a result above the MDL and/or RL, the result has been "B" flagged.

8011 GAS CHROMATOGRAPHY

No difficulties were encountered during the EDB analyses.

SEMIVOLATILE ORGANIC COMPOUNDS (GC-MS)

No difficulties were encountered during the semivolatiles analyses.

GAS RANGE ORGANICS

GRO (C6-C10) was detected in method blank MB 480-30437/3 at a level that was above the method detection limit but below the reporting limit. The value should be considered estimate, and has been flagged "J". If the associated sample reported a result above the MDL and/or RL, the result has been "B" flagged.

DISSOLVED GASES

Due to the concentration of target compounds detected, samples H-1 (180-3644-1)[1000X], FH-1 (180-3644-2)[20X], R-1 (180-3644-5)[1000X], R-1 (180-3644-5)[200X], R-2 (180-3644-6)[100X], R-2 (180-3644-6)[1000X] and KDE-1 (180-3644-9)[10X] were analyzed at a dilution. The reporting limits have been adjusted accordingly.

GLYCOLS

The matrix spike and matrix spike duplicate of sample H-1 (180-3644-1) recovered outside of the control limits for Triethylene Glycol.

DIESEL RANGE ORGANICS

No difficulties were encountered during the DRO analyses.

METALS

The serial dilution of sample H-1 (180-3644-1) was outside of the percent difference control limits for several metals.

The method blanks had analytes detected at levels that were above the method detection limit but below the reporting limit. The values should be considered estimates, and have been flagged "J". If the associated sample reported a result above the MDL and/or RL, the result has been "B" flagged.

The sodium concentrations found in the continuing calibration blanks three thru eight (CCB3 thru CCB8) were greater than the reporting limit (100ug/l). All associated samples bracketed by these CCB's had sodium concentrations at least 100x greater than the CCB concentrations. The CCB concentrations were more than likely caused by "memory effect" and no positive bias in the results is suspected.

The manganese concentrations found in the continuing calibration blanks five, seven and eight (CCB5, CCB7 and CCB8) were greater than the reporting limit (0.5ug/l). All associated samples bracketed by these CCB's had manganese concentrations at least 10x greater than the CCB concentrations. The CCB concentrations were more than likely caused by "memory effect" and no positive bias in the results is suspected.

GENERAL CHEMISTRY

pH is a field parameter. Laboratory pH analysis was completed at the request of the client.

The samples were analyzed outside of the holding time for MBAS.

Due to the matrix, the initial volume used for the following sample deviated from the standard procedure for method 2540D: S-1 (180-3644-4). The reporting limits (RLs) have been adjusted proportionately.

The method blanks had analytes detected at levels that were above the method detection limit but below the reporting limit. The values should be considered estimates, and have been flagged "J". If the associated sample reported a result above the MDL and/or RL, the result has been "B" flagged.