



3058 Research Drive
State College, Pennsylvania 16801 USA
Telephone: 814.272.1039
Fax: 814.272.1019

Analytical Report

PFOA and PFOS Analysis of Deer Muscle and Liver Samples by LC/MS/MS

MPI Report No. L0019346

Testing Laboratory

MPI Research, Inc.
3058 Research Drive
State College, PA 16801

Requester/Project Manager

Dena Haverland
Dalton Utilities
PO BOX 869
Dalton, GA 30722
Phone: 706-529-1010

2010 JMW - 6 P 1:44

1 Introduction

Results are reported for the analysis of deer muscle and liver samples received at MPI Research from Dalton Utilities. The MPI Research study number assigned to the project is L0019346. Table I lists the target analytes quantitated for the samples.

Table I. Target Analytes for Quantitation

Compound Name	Acronym
Perfluorooctanoic Acid	C8 Acid or PFOA
Perfluorooctanesulfonate	C8 Sulfonate or PFOS

2 Sample Receipt

Four samples were received from Dena Haverland at Dalton Utilities for this study. The samples were collected on October 02, 2009. The samples arrived on October 06, 2009 via Fedex and were logged in under MPI Research login number L0019346. The shipment was received frozen on dry ice. The samples were stored frozen at approximately -20°C from receipt until analysis. Chain-of-custody information is presented in Attachment A.

3 Methods - Analytical and Preparatory

3.1 Muscle and Liver Sample Preparation

- 3.1.1. Weigh 1 g of muscle or liver sample into a 50 mL disposable centrifuge tube and fortify, if appropriate. Add 0.2 mL of a 100 ng/mL WIS for a final concentration of 1 ng/mL.
- 3.1.2. Add water to the sample for a final volume of 10 mL. Cap tightly.
- 3.1.3. Homogenize sample using a tissuemizer for ~1 minute.
- 3.1.4. Transfer 1 mL of the sample using a disposable pipette into 15 mL disposable centrifuge tubes. Add 5 mL of ACN and shake for ~20 minutes on a wrist action shaker.
- 3.1.5. Centrifuge tubes at ~3000 rpm for ~ 5 minutes. Carefully decant supernatant into a 50 mL disposable centrifuge tube and add 35 mL of water.
- 3.1.6. Place the unconditioned SPE columns on the vacuum manifold. Condition the SPE columns by passing ~ 10 mL of methanol through the column followed by ~ 5 mL of water. The washes may be pulled through the SPE column using vacuum at a flow rate of ~1 drop/sec or may be allowed to pass through the column unaided. Discard all washes. Do not allow the column to dry.
- 3.1.7. Load the sample onto a conditioned SPE column . Discard the eluate. Any analyte residues will be trapped on the SPE column at this point.
- 3.1.8. Elute with 2 mL of methanol. Collect 2 mL of elute into a graduated 15 mL centrifuge tube.

3.2 Sample Analysis by LC/MS/MS

In High Pressure Liquid Chromatography (HPLC), an aliquot of extract is injected and passed through a liquid-phase chromatographic column. Based on the affinity of the analyte for the stationary phase in the column relative to the liquid mobile phase, the analyte is retained for a characteristic amount of time. Following HPLC separation, mass spectrometry provides a rapid and accurate means for analyzing a wide range of organic compounds. Molecules are ionized, fragmented, and detected. The ions characteristic of the compounds are observed and quantitated against external calibration standards.

An HP1100 system interfaced to an Applied Biosystems API 4000 LC/MS/MS was used to analyze the sample extracts for quantitation. A gradient elution through a Phenomenex Luna 3 μ C8(2) Mercury, 20 x 4.0 mm column was used for separation.

The following gradient was performed:

Mobile Phase (A): 2mM Ammonium Acetate in Water
Mobile Phase (B): Methanol

Time	%A	%B
0.0	90	10
0.5	90	10
2.0	10	90
5.0	10	90
5.1	0	100
6.0	0	100
6.1	90	10
10.0	90	10

The following parameters were used for operation of the mass spectrometer:

Parameter	Setting
Ionization Mode	Electrospray
Polarity	Negative
Transitions Monitored	413→369 (PFOA) 499→80 (PFOS) 415→370 (Internal Std. ^{13}C PFOA (m+2)) 503→80 (Internal Std. ^{13}C PFOS (m+4))
Gas Temperature	450°C

4 Analysis by LCMSMS

4.1 Calibration

For the muscle and liver sample analysis, a 6-point calibration curve was analyzed throughout the analytical sequence for PFOA and PFOS. The calibration points were prepared at 0.1, 0.2, 0.5, 1.0, 2.0, 5.0 ng/mL (ppb) each containing 1.0 ng/mL ^{13}C -PFOA (m+2) and ^{13}C -PFOS (m+4).

The ratio of the analyte concentration to the IS concentration versus the ratio of the analyte instrument response (area) to the IS response (area) was plotted for each point. Using linear regression with 1/x weighting, the slope, y-intercept and coefficient of determination (r^2) were determined. A calibration curve is acceptable if $r^2 \geq 0.985$.

For the results reported here, calibration criteria were met. The calibration curves are included in the raw data in Attachment C.

4.2 Laboratory Control Spikes

Laboratory control spikes in the analytical set were prepared during each extraction set by adding a known concentration of the analyte to deer muscle and liver controls. Laboratory control spikes are used to assess method accuracy. The laboratory control spikes must show recoveries between 70-130% or the data is rejected. For the results reported here, the laboratory control spikes were within the acceptable range. Laboratory control spike recoveries are given in Attachment B.

4.3 Matrix Spikes

A matrix spike was prepared for each sample by adding a known concentration of the target analyte to a sample. Matrix spikes are used to assess method accuracy in the matrix. The matrix spikes should show recoveries between 70-130%. For the results reported here, the matrix spike was within the acceptable range with the exceptions of:

L19346-4 (Deer # 7 3.5 yr male liver) Spk C at 2000 ng/g for PFOS, which gave low recovery of 59%.

4.4 Laboratory Duplicates

Each sample was prepared in duplicate and analyzed. Duplicate results are given along with the sample results in Attachment B.

5 Data Summary

Please see Attachment B for a detailed listing of the analytical results. For the muscle and liver samples the results are reported in parts per billion (ng/g) on an as-received basis.


6 Data/Sample Retention

Samples are disposed of 60 days after the report is issued unless otherwise specified by the project manager. All electronic data is archived on retrievable media and hard copy reports are stored in data folders maintained by MPI Research. Hardcopy data is stored for a minimum of five years. The client will be notified 30 days prior to the disposal of hardcopy data.

7 Attachments

- 7.1 Attachment A: Chain of Custody
- 7.2 Attachment B: Analytical Results
- 7.3 Attachment C: Raw Analytical Data for Water


8 Signatures



Mark Neeley, Research Chemist Associate II

11-6-09

Date



Robert Zhu, Manager, Analytical

11/8/05

Date

A

Login

Login Group: L0019346

Login #:	19460	Conform COC Sample:	True
Project:	P0005196	Conform COC:	True
Company Name:	Dalton Utilities	Conform Sample:	True
Submitted By:	Dena Haverland	Conform Request:	True
Login Type:	Immediate Receipt of Samples		
Started:	True		
Date Start:	10/27/2009		
Due Date:	11/06/2009		
Login Initiated:	10/27/2009		
Received By:	Ammerman, Mark		
Spread Sample:			
Label:			
MPI SD/PI:	Zhu, Xiang		
Project Title/Type: PFOA and PFOS Analysis of Animal Muscle and Liver by LC/MS/MS / ROUTINE			
Login Notes:			

Packages / Containers

<u>Package</u>	<u>Carton</u>	<u>Date / Condition</u>		<u>Shipper / ID</u>	<u>Temp. Control/Temp.</u>	<u>Direction / Handled By</u>
K0022042		Received Date: 10/27/09 10:25 Package & Contents Uncompromised		FEDEX 8694 2057 8178	Dry Ice -79.2	RECEIVED Ammerman, Mark
Container #	Gross Weight	pH	Container Type	Preservative	Mfg. Lot	Mfg. ID
C0457624	218.20 g		1/2 gallon ziplock bag	NONE		
C0457625	278.00 g		1/2 gallon ziplock bag	NONE		
C0457626	324.10 g		1 gallon ziploc bag	NONE		
C0457627	994.10 g		1 gallon ziploc bag	NONE		

Samples

<u>Sample ID</u>	<u>Container</u>	<u>Matrix</u>	<u>System</u>	<u>System Matrix</u>	<u>Sample</u>	<u>Date Sampled</u>	<u>Date Due</u>
L0019346-0001	C0457624	SOLID	Deer	Tissue	Deer #6 0.5 yr female-muscle	10/02/2009	11/06/2009
L0019346-0002	C0457626	SOLID	Deer	Liver	Deer #6 0.5 yr female-liver	10/02/2009	11/06/2009
L0019346-0003	C0457625	SOLID	Deer	Tissue	Deer #7 3.5 yr male-muscle	10/02/2009	11/06/2009
L0019346-0004	C0457627	SOLID	Deer	Liver	Deer #7 3.5 yr male-liver	10/02/2009	11/06/2009



Login

Login Reviewed By:

RHCL

Date/Time:

12/28/09 1440



Sample Submittal



MPI Research Contact: Daniel Wright

Send Report To:	
Company:	Dalton Utilities 1200 VD Parrott JR Parkway, PO Box 869
Address:	
City, State, ZIP:	Dalton, GA 30722-0869
Attention:	Dena Haverland
Phone #:	706-529-1010
Fax #:	706-529-1271
Email:	dahaverland@dutil.com
Study/Job #:	
Signature/Date:	
Printed Name:	

Please fax this form before sending samples.

Please send samples to shipping and receiving:
3048 Research Drive, State College, PA 16801
T: (814) 272-1039 • F: (814) 272-1019

Turnaround time (TAT) requirements:

Results Due Date: 30 days
 Preliminary Results Format: Verbal Email Fax
 Report Due Date: 30 days

Storage Conditions	Safety Information
Room temperature Refrigerator <input checked="" type="checkbox"/> Freezer <input checked="" type="checkbox"/> Ultra Low freezer Desiccated Lighting required Stability (%C/%RH): Stability time period:	Special handling: _____ MSDS attached Controlled substance: _____ HAZARDS: _____ Please fill in the diamond HMIS/NFPA (0-4) if appropriate

Client ID# Description	Lot/ Control #	Amt. Sent/ Weight	# of Bottles	Matrix	Date & Time	Tests Requested
1 Deer #6 0.5 yr female - Serum		10 ml	10	deer	10/2/09 1:08AM	PFOA/PFOS
2 Deer #6 0.5 yr female - muscle		requested	1 bag	deer	10/2/09 2:28AM	PFCA/PFOS
3 Deer #6 0.5 yr female - Liver		whole	1 bag	deer	10/2/09 2:36AM	PFOA/PFOS
4 Deer #7 3.5 yr Male - Serum		10 ml	10	deer	10/2/09 1:45AM	PFOA/PFOS
5 Deer #7 3.5 yr Male - Muscle		requested	1 bag	deer	10/2/09 2:45AM	PFCA/PFOS
6 Deer #7 3.5 yr Male - Liver		whole	1 bag	deer	10/2/09 2:48AM	PFOA/PFOS
7						
8						
9						
10						

PO #

Notes

Relinquished by: Date: Time: Received by: Date: Time:
 Daniel Wright 10/5/09 6:35AM 10/6/09 10:27

'THIS IS AN EXACT COPY OF
THE ORIGINAL DOCUMENT'

I NY MW DATE 10/27/09



TEMPORARY SAMPLE STORAGE FORM

To be completed during ExyLIMS Login

Project #: P 5196

Login #: L 4346

Initials / Date: MCA 10/27/09

One form to be completed for each package

Date / Time Received: 10/26/09 1025

Received By: Mark Amerson

Shipper: FedEx

Shipper Package ID: 8694 2057 8178

Temperature (deg C) / Thermometer ID: -79.2 / D0001775

Temperature Control Method: dry ice/ice

Temporary Storage Location: walk-in freezer 11

Condition of sample(s):

- Good – Package and contents uncompromised
- Fair – Package damaged / contents uncompromised
- Poor – Package and contents compromised

Notes:

FedEx® USAirbill

FedEx
Tracking
Number

8694 2057 8178



1 From	2 To	3 To	4a Express Package Service	4b Express Freight Service
From <i>10/5/01</i>	To <i>300 W. 1st Street</i>	Recipient's Name <i>John Doe</i>	<input checked="" type="checkbox"/> FedEx Priority Overnight <small>Next business day delivery. Not available Saturday. Delivery is restricted.</small>	<input type="checkbox"/> FedEx Standard Overnight <small>Next business day delivery.</small>
Date <i>10/5/01</i>	Phone <i>(319) 344-1234</i>	Name <i>John Doe</i>	<input type="checkbox"/> FedEx 2 Day/ <small>Second business day* Thursday shipments will be delivered on Monday unless Saturday Delivery is selected.</small>	<input type="checkbox"/> FedEx First Overnight <small>Same day delivery. Not available Saturday. Delivery NOT available.</small>
Sender's Name <i>John Doe</i>	Address <i>123 Main Street</i>	Company <i>ABC Corp.</i>	<input type="checkbox"/> FedEx Express Saver <small>Third business day* Saturday Delivery NOT available.</small>	<input type="checkbox"/> FedEx 3 Day/Freight <small>Second business day* Thursday shipments will be delivered on Monday unless Saturday Delivery is selected.</small>
City <i>Iowa City</i>	State <i>IA</i>	ZIP <i>31930</i>	<small>* FedEx Overnight rate not available. Minimum charge One-pound rate.</small>	<small>* To most locations.</small>
Dest. Phone/Store <i>(319) 344-1234</i>			Packages up to 150 lbs.	

2 Your Internal Billing Reference

City <i>Iowa City</i>	State <i>IA</i>	ZIP <i>31930</i>
Dest. Phone/Store <i>(319) 344-1234</i>		
Date/Pickup/Station <i>10/5/01</i>		

Address

To receive a package held at a specific FedEx location, print FedEx address here.



8694 2057 8178

3 To Your Internal Billing Reference

Recipient's Name <i>John Doe</i>	Phone <i>(319) 344-1234</i>
Address <i>300 W. 1st Street</i>	
Company <i>ABC Corp.</i>	
Recipient's Address <i>300 W. 1st Street</i>	
We cannot deliver to P.O. boxes or P.O. ZIP codes.	
Address	

5 Packaging

- FedEx Envelope* FedEx Pak*
Includes FedEx Seal Pak, FedEx Super Pak, FedEx Lumper Pak, and FedEx Sure Pak.

6 Special Handling

- SATURDAY Delivery Not available for FedEx Same Day, Friday Night, Saturday, FedEx Express, FedEx Ground, FedEx Home Delivery, FedEx Smart Cut, FedEx Direct, FedEx Next Day, FedEx Next Day Overnight.
- HOLD Location Not available for FedEx Same Day, Friday Night, Saturday, FedEx Express, FedEx Ground, FedEx Home Delivery, FedEx Smart Cut, FedEx Direct, FedEx Next Day, FedEx Next Day Overnight.
- HOLD Saturday Available ONLY for FedEx Priority Overnight and FedEx Two Day.
- HOLD Sunday Available ONLY for FedEx Priority Overnight and FedEx Two Day.
- No One hour must be checked.
- Yes As per attached Shipper's Declaration, not required.
- Dangerous goods including dry ice cannot be shipped in FedEx products.

7 Payment/Bill to:

- Sender Enter FedEx Acct. No. or Credit Card No. below.
- Recipient Or obtain Recip. Acct. No. Credit Card Cash/Check
- Third Party Credit Card Cash/Check

Total Packages Total Weight Total Declared Value

Your liability is limited to \$100 unless you declare a higher value. See back for details.

Credit Card Auth.

8 Residential Delivery Signature Options

If you require a signature, check Direct or Indirect.

- Required Direct Signature Indirect Signature
- FedEx may be held Someone at recipient's address may sign for delivery. New address may apply.
- FedEx may be held Indirect Signature

520

fedex.com 1800 GoFedEx 1800.463.3393

IS AN EXACT COPY OF
ORIGINAL DOCUMENT

BY

M/A DATE *10/27/01*

B



3058 Research Drive
State College, Pennsylvania 16801 USA
Telephone: 814.272.1039
Fax: 814.272.1019

Analytical Report

Summary of Fluorochemical Residues in Muscle Samples

Sample ID	PFOA Perfluorooctanoic Acid	PFOS Perfluorooctanesulfonate
	Analyte Found (ng/g, ppb)	Analyte Found (ng/g, ppb)
Deer # 7 3.5 yr male-muscle	ND	88.4
Deer # 7 3.5 yr male-muscle*	ND	92.2
Deer # 6 0.5 yr female-muscle	ND	13.4
Deer # 6 0.5 yr female-muscle*	ND	11.9

*Laboratory Duplicate

ND = Not detected = Response is below the LOD of 1.0 ng/g (ppb).

NQ = Not quantifiable = Response is between the LOD and the LOQ of 10 ng/g (ppb).



3058 Research Drive
State College, Pennsylvania 16801 USA
Telephone: 814.272.1039
Fax: 814.272.1019

Analytical Report

Recovery Summary of Fluorochemical Residues in Muscle Samples

Sample Description	Amount Spiked (ng/g)	PFOA			PFOS		
		Amt Found in Sample (ng/g)	Amount Recovered (ng/g)	Recovery (%)	Amt Found in Sample (ng/g)	Amount Recovered (ng/g)	Recovery (%)
LCS A (Data set 110409A) 10 ng/g	10	ND	9.37	94	ND	8.78	88
LCS B (Data set 110409A) 50 ng/g	50	ND	47.1	94	ND	45.2	90
Deer # 7 3.5 yr male-muscle (L19346-3 Spk C, 50 ng/g Lab Spike)	50	ND	48.5	97	88	128	79
Deer # 6 0.5 yr female-muscle (L19346-1 Spk D, 50 ng/g Lab Spike)	50	ND	46.4	93	13.4	58.0	89

ND = Not detected = Response is below the LOD of 1.0 ng/g.

NQ = Not quantifiable = Response is between the LOD and the LOQ of 10 ng/g.



3058 Research Drive
State College, Pennsylvania 16801 USA
Telephone: 814.272.1039
Fax: 814.272.1019

Summary of Fluorochemical Residues in Liver Samples

Sample ID	PFOA Perfluorooctanoic Acid	PFOS Perfluorooctanesulfonate
	Analyte Found (ng/g, ppb)	Analyte Found (ng/g, ppb)
Deer # 7 3.5 yr male-liver	ND	2150^
Deer # 7 3.5 yr male-liver*	ND	2160^
Deer # 6 0.5 yr female-liver	ND	604
Deer # 6 0.5 yr female-liver*	ND	574

*Laboratory Duplicate

ND = Not detected = Response is below the LOD of 1.0 ng/g (ppb).

NQ = Not quantifiable = Response is between the LOD and the LOQ of 10 ng/g (ppb).

[^] The corresponding Matrix Spike recovery was outside the acceptance criteria of 70-130%, therefore the sample values should be considered an estimate.



3058 Research Drive
State College, Pennsylvania 16801 USA
Telephone: 814.272.1039
Fax: 814.272.1019

Recovery Summary of Fluorochemical Residues in Liver Samples

Sample Description	Amount Spiked (ng/g)	PFOA			PFOS		
		Amt Found In Sample (ng/g)	Amount Recovered (ng/g)	Recovery (%)	Amt Found In Sample (ng/g)	Amount Recovered (ng/g)	Recovery (%)
LCS A (Data set 110409B) 10 ng/g	10	ND	10.0	100	16.5	26.3	98
LCS B (Data set 110409B) 50 ng/g	50	ND	50.3	101	16.5	61.6	90
LCS A (Data set 110509B) 10 ng/g	10	N/A	N/A	N/A	18.6	28.8	102
LCS B (Data set 110509B) 50 ng/g	50	N/A	N/A	N/A	18.6	67.9	99
Deer # 7 3.5 yr male-liver (L19346-4 Spk C, 50 ng/g Lab Spike)	50	ND	52.4	105	2150	**	**
Deer # 6 3.5 yr female-liver (L19346-2 Spk D, 50 ng/g Lab Spike)	50	ND	50.4	101	604	**	**
Deer # 7 3.5 yr male-liver (L19346-4 Spk C, 2000 ng/g Lab Spike)	2000	N/A	N/A	N/A	2150	3330	59*
Deer # 6 3.5 yr female-liver (L19346-2 Spk D, 250 ng/g Lab Spike)	250	N/A	N/A	N/A	604	817	85

ND = Not detected = Response is below the LOD of 1.0 ng/g.

NQ = Not quantifiable = Response is between the LOD and the LOQ of 10 ng/g.

* The Matrix Spike recovery was outside the acceptance criteria of 70-130%, therefore the sample values should be considered an estimate.

** The endogenous level of PFOS in the sample significantly exceeds the spiking level, therefore an accurate recovery cannot be calculated.

C

RAW DATA REPORT

Sponsor Study No: NA Limit of Quantitation: 10 ng/g Set No: 110409A
 MPI Study No: L19345 Injection Volume: 15 µL Analyst: Mark Neeley
 Analyte: PFOA Matrix: Deer Muscle Instrument Type: LC/MS/MS Unit # 9
 Ions Monitored: 413 -> 369 Sample Weight 1.0 g Extraction Date: 11/04/09
 Site: NA Analyzed on: 11/04/09

MPI Research ID	Sponsor ID	Sample Code	Sample Index No.	Internal					Internal		Amount	
				Std. Conc. (ng/mL)	Std. Conc. (ng/mL)	Aliquot Factor (AF)	Dilution Factor (DF)	Peak Area	Standard Peak Area	Analyte Found (ng/g)	Analyte Added (ng/g)	Recovery (%)
SS33618	-	CS	1	0.100	1.0	-	-	55054	590483	-	-	-
SS33617	-	CS	2	0.200	1.0	-	-	104150	572151	-	-	-
SS33616	-	CS	3	0.500	1.0	-	-	253254	591894	-	-	-
SS33615	-	CS	4	1.00	1.0	-	-	512706	573466	-	-	-
SS33614	-	CS	5	2.00	1.0	-	-	1006687	578731	-	-	-
SS33613	-	CS	6	5.00	1.0	-	-	2362116	567936	-	-	-
Methanol Wash	-	W	7	-	-	-	-	1476	0	-	-	-
Methanol Wash	-	W	8	-	-	-	-	1517	0	-	-	-
Control	MC4314 Deer Muscle Control	C	9	-	1.0	20	1	12498	630414	ND	-	-
LCS A	MC4314 Deer Muscle Spike A	LCS	10	-	1.0	20	1	311604	765110	9.37	10	94
LCS B	MC4314 Deer Muscle Spike B	LCS	11	-	1.0	20	1	1418119	709162	47.1	50	94
L19346-3 Spk C	Deer # 7 3.5 yr male-muscle Spike C	LF	12	-	1.0	20	1	1437886	698422	48.5	50	97
L19346-1 Spk D	Deer # 6 0.5 yr female-muscle Spike D	LF	13	-	1.0	20	1	862089	438366	46.4	50	93
L19346-3	Deer # 7 3.5 yr male-muscle	S	14	-	1.0	20	1	5207	670093	ND	-	-
L19346-3 Dup	Deer # 7 3.5 yr male-muscle Duplicate	S	15	-	1.0	20	1	5129	649587	ND	-	-
L19346-1	Deer # 6 0.5 yr female-muscle	S	16	-	1.0	20	1	3752	506345	ND	-	-
L19346-1 Dup	Deer # 6 0.5 yr female-muscle Duplicate	S	17	-	1.0	20	1	3129	443379	ND	-	-
SS33615		CCV	18	1.00	1.0	-	-	525395	590143	1.04	1.0	104

Analyte Found (ng/g) = (((analyte peak area/IS peak area) - intercept) / slope) x IS conc.) x AF x DF/Sample weight

Standard Curve: Linear (1/x weighted)

Intercept = 0.0125

Slope = 0.843

Coef. Of Det. = 0.9992

Recovery (%) = $\frac{[\text{Analyte found (ng/g)} - \text{Analyte found in control (ng/g)}]}{\text{amount Analyte added (ng/g)}} \times 100$

CS = Calibration standard

LF = Lab fortified sample

W = Methanol Wash

CCV = Continuing Calibration Verification

FF = Field fortified sample

ND = Not detected = Response between 0 and LOD

C = Control sample

LCS = Laboratory Control Spike

NQ = Not quantifiable = Response between LOD and LOQ

Spreadsheet prepared by: Mark / 11-6-09

RAW DATA REPORT

Sponsor Study No:	NA	Limit of Quantitation:	10 ng/g	Set No:	110409A
MPI Study No:	L19345	Injection Volume:	15 µL	Analyst:	Mark Neetley
Analyte:	PFOS	Matrix:	Deer Muscle	Instrument Type:	LC/MS/MS Unit # 9
Ions Monitored:	499 > 80	Sample Weight	1.0 g	Extraction Date:	11/04/09
Site:	NA			Analyzed on:	11/04/09

MPI Research ID	Sponsor ID	Sample Code	Sample Index No.	Internal				Internal		Amount		
				Std. Conc. (ng/mL)	Std. Conc. (ng/mL)	Extraction Factor	Dilution Factor	Peak Area	Peak Area	Standard Found (ng/g)	Analyte Added (ng/g)	Recovery (%)
SS33618	-	CS	1	0.100	1	-	-	39149	344167	-	-	-
SS33617	-	CS	2	0.200	1	-	-	75056	338773	-	-	-
SS33616	-	CS	3	0.500	1	-	-	187236	342440	-	-	-
SS33615	-	CS	4	1.00	1	-	-	369109	342687	-	-	-
SS33614	-	CS	5	2.00	1	-	-	747914	338050	-	-	-
SS33613	-	CS	6	5.00	1	-	-	1751107	332437	-	-	-
Methanol Wash	-	W	7	-	-	-	-	0	0	-	-	-
Methanol Wash	-	W	8	-	-	-	-	0	0	-	-	-
Control	MC4314 Deer Muscle Control	C	9	-	1	20	1	2080	301876	ND	-	-
LCS A	MC4314 Deer Muscle Spike A	LCS	10	-	1	20	1	182837	381538	8.78	10	88
LCS B	MC4314 Deer Muscle Spike B	LCS	11	-	1	20	1	852146	350892	45.2	50	90
L19346-3 Spk C	Deer # 7 3.5 yr male-muscle Spike C	LF	12	-	1	20	1	2466358	330184	*	50	*
L19346-1 Spk D	Deer # 6 0.5 yr female-muscle Spike D	LF	13	-	1	20	1	710987	228483	58.0	50	89
L19346-3	Deer # 7 3.5 yr male-muscle	S	14	-	1	20	1	1572929	332029	88.4	-	-
L19346-3 Dup	Deer # 7 3.5 yr male-muscle Duplicate	S	15	-	1	20	1	1534874	310563	92.2	-	-
L19346-1	Deer # 6 0.5 yr female-muscle	S	16	-	1	20	1	183224	251539	13.4	-	-
L19346-1 Dup	Deer # 6 0.5 yr female-muscle Duplicate	S	17	-	1	20	1	153884	237149	11.9	-	-
SS33615	-	CCV	18	1.00	1	-	-	369522	346418	0.988	1.00	99

Analyte Found (ppb) = (((analyte peak area/IS peak area) - intercept) / slope) x IS conc. x DF

Standard Curve: Linear (1/x weighted)

Intercept = 0.00957

Slope = 1.07

Coef. Of Det. = 0.9996

Recovery (%) = $\frac{[\text{Analyte found (ng/g)} - \text{Analyte found in control (ng/g)}]}{\text{amount Analyte added (ng/g)}} \times 100$

CS = Calibration standard

LF = Lab fortified sample

W = Methanol Wash

CCV = Continuing Calibration Verification

FF = Field fortified sample

ND = Not detected = Response between 0 and LOD

C = Control sample

LCS = Laboratory Control Spike

NQ = Not quantifiable = Response between LOD and LOQ

* Sample response exceeds the calibration range, dilution required, see set no 110409AD for analysis

Spreadsheet prepared by: MDN/ 11-6-09



3058 Research Drive
State College, PA 16801

Phone: 814-272-1039
Fax: 814-231-1580

Internal Chain of Custody/Fortification Sheet

MPI Research Study Number:
Sponsor Study/Protocol No:

L19346
NA

Matrix: Deer Muscle

The samples listed below were removed from Freezer No. 36
Time 9:25 AM

Date 11-4-09

Initials MAR

CLIENT SAMPLE ID	MPI RESEARCH ID NUMBER	WEIGHT (g)	FORTIFICATION (ng)
na	MC4314 Control	1.0	-
na	MC4314 LCS A	1.0	10
na	MC4314 LCS B	1.0	50
Deer #7 3.5 yr male-muscle	L19346-3 Spk C	1.0	50
Deer #6 0.5 yr female-muscle	L19346-1 Spk D	1.0	50
Deer #6 0.5 yr female-muscle	L19346-1	1.0	-
Deer #6 0.5 yr female-muscle	L19346-1 Dup	1.0	-
Deer #7 3.5 yr male-muscle	L19346-3	1.0	-
Deer #7 3.5 yr male-muscle	L19346-3 Dup	1.0	-
-	-	-	-
-	-	-	-
-	-	-	-
-	-	-	-
-	-	-	-
-	-	-	-
-	-	-	-
-	-	-	-
-	-	-	-
-	-	-	-
-	-	-	-
-	-	-	-
-	-	-	-

	Spiking Solution Used	Volume Used for Spiking	Initial/Date
MC4314 LCS A	SS0033620 (200 ng/mL)	50 µL (50-250 µL auto-pipette)	<u>MAR</u> / <u>11-4-09</u>
MC4314 LCS B	SS0033619 (1000 ng/mL)	50 µL (50-250 µL auto-pipette)	<u>MAR</u> / <u>11-4-09</u>
L19346-3 Spk C	SS0033619 (1000 ng/mL)	50 µL (50-250 µL auto-pipette)	<u>MAR</u> / <u>11-4-09</u>
L19346-1 Spk D	SS0033619 (1000 ng/mL)	50 µL (50-250 µL auto-pipette)	<u>MAR</u> / <u>11-4-09</u>
All Samples	SS0032652 (100 ng/mL)	200 µL (50-250 µL auto-pipette)	<u>MAR</u> / <u>11-4-09</u>

All samples were weighed on balance 20.

Time 9:26 AM

Date 11-4-09

Initials MAR

After weighing samples were returned to Freezer No. 36

Time 9:39 AM

Date 11-4-09

Initials MAR

Comments:

Analysis Summary:

Data Set: 110409A
Data Set: 110409AD
Data Set: -

Initials/Date: MAR / 11-4-09
Initials/Date: MAR / 11-5-09
Initials/Date: - / -

Set extraction/analysis data verified by: Amb

Date: 11/09/09

Nov 02, 2009/2

Object: \\sc1wp5556\mdrive\PE SCIEX DATA\Projects\P5195 Batch:09_110409A Muscle Tab:Sample Set:SET1 AcqMethod:P5195_102
m/n 11-479

Sam

Sample Name	Sample ID	Vial Position	Data File
SS33618	Calibration Standard, 0.1 ng/mL	1	09_110409A\110409A
SS33617	Calibration Standard, 0.2 ng/mL	2	09_110409A\110409A
SS33616	Calibration Standard, 0.5 ng/mL	3	09_110409A\110409A
SS33615	Calibration Standard, 1.0 ng/mL	4	09_110409A\110409A
SS33614	Calibration Standard, 2.0 ng/mL	5	09_110409A\110409A
SS33613	Calibration Standard, 5.0 ng/mL	6	09_110409A\110409A
Methanol Wash	Methanol Wash	91	09_110409A\110409A
Methanol Wash	Methanol Wash	91	09_110409A\110409A
Control	MC4314 Deer Muscle Control	51	09_110409A\110409A
LCS A	MC4314 Deer Muscle Spike A, 10 ng/mL	52	09_110409A\110409A
LCS B	MC4314 Deer Muscle Spike B, 50 ng/mL	53	09_110409A\110409A
L19346-3 Spk C	Deer # 7 3.5 yr male-muscle Spike C, 50 ng/mL	54	09_110409A\110409A
L19346-1 Spk D	Deer # 6 0.5 yr female-muscle Spike D, 50 ng/mL	55	09_110409A\110409A
L19346-3	Deer # 7 3.5 yr male-muscle	56	09_110409A\110409A
L19346-3 Dup	Deer # 7 3.5 yr male-muscle Duplicate	57	09_110409A\110409A
L19346-1	Deer # 6 0.5 yr female-muscle	58	09_110409A\110409A
L19346-1 Dup	Deer # 6 0.5 yr female-muscle Duplicate	59	09_110409A\110409A
SS33615	CCV, 1.0 ng/mL	4	09_110409A\110409A

LC/MS/MS SYSTEM AND OPERATING CONDITIONS

Protocol No: NA

MPI Study No: L19346

Instrument: AB API 4000 Biomolecular Mass Analyzer, (LC/MS/MS #9)
SCIEX Turbo Ion Spray Liquid Introduction Interface
Turbo Ion spray temperature = 450 °C

Computer: Dell OptiPlex GX 110

Software: PE Sciex Analyst 1.4

HPLC Equipment: Hewlett Packard (HP) Series 1100
HP Quat Pump HP Vacuum Degasser
HP Autosampler HP Column Oven

HPLC Column: Phenomenex Luna C8 (2) Mercury, 2cm x 4mm, 3 µm (ExyLIMS ID:
MA0052622)

Column Temperature: 35°C

Mobile Phase (A): 2 mM Ammonium Acetate in Water (ExyLIMS ID: SL0045925)

Mobile Phase (B): Methanol (ExyLIMS ID: RE0047880)

Injected Volume: 15 µL

Time (min)	% A	% B	Flow Rate (µL/min)
0.0	90	10	750
0.5	90	10	750
2.0	10	90	750
5.0	10	90	750
5.1	0	100	750
6.0	0	100	750
6.1	90	10	750
10.0	90	10	750

Ions monitored:

Analyte	Parent ion	Daughter ion(s)	Dwell (secs)
PFOA	413	369	0.200
PFOS	499	80	0.200
¹³ C PFOA (m+2)	415	370	0.200
Internal Standard			
¹³ C PFOS (m+4)	503	80	0.200
Internal Standard			

Analyst: Mark Neeley *MNP 11-4-09*
MPI Research, Inc.
3058 Research Drive, State College, PA 16801
Phone: (814) 272-1039 FAX: (814) 231-1580

All Handwritten Peak ID's by: *MNP 11-5-09*

alyst Version: 1.4.2
rinting Time: 1:50:07 PM
rinting Date: Wednesday, November 04, 2009

MPI Study: L19346
MPI Set No.: 110409A

Operator: Mark Neeley
Instrument No.: LC/MS/MS #9

MDP 11-4-09

Acquisition Information:

quisition Method: P5195_102909.dam
reated: Thursday October 29 2009 09: 35: 12 AM
Last Modified: Friday October 30 2009 10: 56: 20 AM
Comment: PFOA PFOS
ynchronization Mode: LC Sync
uto-Equilibration: Off
Acquisition Duration: 10min0sec
Number Of Scans: 732
eriods In File: 1
Acquisition Module: Acquisition Method
Software version Analyst 1.4.2

Period 1:

Scans in Period: 732
relative Start Time: 0.00 msec
Experiments in Period: 1

Period 1 Experiment 1:

Scan Type: MRM (MRM)
Polarity: Negative
Scan Mode: N/A
Ion Source: Turbo Spray
Resolution Q1: Unit
Resolution Q3: Unit
Intensity Thres.: 0.00 cps
Setting Time: 0.0000 msec
MR Pause: 5.0070 msec
ICA: No
Step Size: 0.00 amu

Q1 Mass (amu)	Q3 Mass (amu)	Dwell (msec)	Param	Start	Stop
13.00	369.00	200.00	DF	-32.00	-32.00
			CE	-18.00	-18.00
			CXF	-13.20	-13.20

Q1 Mass (amu)	Q3 Mass (amu)	Dwell (msec)	Param	Start	Stop
415.00	370.00	200.00	DF	-32.00	-32.00
			CE	-18.00	-18.00
			CXF	-13.20	-13.20

Q1 Mass (amu)	Q3 Mass (amu)	Dwell (msec)	Param	Start	Stop
499.00	80.00	200.00	DF	-83.00	-83.00
			CE	-88.00	-88.00
			CXF	-6.00	-6.00

Q1 Mass (amu)	Q3 Mass (amu)	Dwell (msec)	Param	Start	Stop
503.00	80.00	200.00	DF	-83.00	-83.00
			CE	-88.00	-88.00

alyst Version: 1.4.2
rinting Time: 1:50:07 PM
rinting Date: Wednesday, November 04, 2009

MPI Study: L19346
MPI Set No.: 110409A

Operator: Mark Neeley
Instrument No.: LC/MS/MS #9

MDR 11-4-09

CXF -6.00 -6.00

parameter Table (Period 1 Experiment 1):

AD: 7.00
CUR: 20.00
GS1: 50.00
S2: 40.00
S: -4500.00
TEM: 450.00
the: ON
F: -10.00

Agilent LC Pump Method Properties

Pump Model: Agilent 1100 LC Quaternary Pump
Minimum Pressure (psi): 0.0
Maximum Pressure (psi): 5801.0
Compressibility: 100.0
Dead Volume (μl): 40.0
Stroke Volume (μl): -1.0
Maximum Flow Ramp (ml/min²): 100.0
Maximum Pressure Ramp (psi/sec): 290.0

Step Table:

Step	Total Time (min)	Flow Rate(μl/min)	A (%)	B (%)	C (%)	D (%)	TE#1	TE#2	TE#3	TE#4
0	0.00	750	0.0	0.0	90.0	10.0	open	open	open	open
	0.50	750	0.0	0.0	90.0	10.0	open	open	open	open
	2.00	750	0.0	0.0	10.0	90.0	open	open	open	open
3	5.00	750	0.0	0.0	10.0	90.0	open	open	open	open
4	5.10	750	0.0	0.0	0.0	100.0	open	open	open	open
	6.00	750	0.0	0.0	0.0	100.0	open	open	open	open
	6.10	750	0.0	0.0	90.0	10.0	open	open	open	open
7	10.00	750	0.0	0.0	90.0	10.0	open	open	open	open

Primary Flow Rate (μl/min): 200.0

Flow Sensor Calibration Table Index: 0

Agilent Column Oven Properties

Left Temperature (°C): 35.00
Right Temperature (°C): 35.00
Temperature Tolerance +/- (°C): 1.00
Start Acquisition Tolerance +/- (°C): 0.50

Time Table (Not Used)

Column Switching Valve Installed

Position for first sample in the batch: Left

Same position for all samples in the batch

Agilent Autosampler Properties

Autosampler Model: Agilent 1100 Wellplate Autosampler
Syringe Size (μl): 100
Injection Volume (μl): 15.00
Draw Speed (μl/min): 200.0
Inject Speed (μl/min): 200.0
Needle Level (mm): 0.00
Temperature Control Enabled
Setpoint (4 - 40 C): 4
Push Location: Flush Port
Push Time (1 - 999 sec): 10

Automatic Delay Volume Reduction Not Used

alyst Version: 1.4.2
rinting Time: 1:50:07 PM

MPI Study: L19346
MPI Set No.: 110409A

Operator: Mark Neeley
Instrument No.: LC/MS/MS #9

rinting Date: Wednesday, November 04, 2009

MDP 11-4-09

Equibration Time (sec): 2
Enz Vial/Well Bottom Sensing No
se Custom Injector Program No


alyst Version: 1.4.2
rinting Time: 8:29:25 AM
rinting Date: Thursday, November 05, 2009

MPI Study: L19346
MPI Set No.: 110409A

Operator: Mark Neeley
Instrument No.: LC/MS/MS #9

MDN 11-5-09

■ 09_110409A.rdb (PFOA): "Linear" Regression ("1 / x" weighting): $y = 0.843 x + 0.0125$ ($r = 0.9996$)



alyst Version: 1.4.2
rinting Time: 8:30:11 AM
rinting Date: Thursday, November 05, 2009


MPI Study: L19346
MPI Set No.: 110409A

Operator: Mark Neeley
Instrument No.: LC/MS/MS #9

Sample Name: "SS33618" Sample ID: "Calibration Standard, 0.1 ng/mL" File: "110409A.wif"

Peak Name: "PFOA" Mass(es): "413.0/369.0 amu"


Comment: "Annotation: "



Sample Name: "SS33618" Sample ID: "Calibration Standard, 0.1 ng/mL" File: "110409A.wif"

Peak Name: "13CPFOA(M+2)(IS)" Mass(es): "415.0/370.0 amu"


Comment: "Annotation: "



Sample Name: "SS33617" Sample ID: "Calibration Standard, 0.2 ng/mL" File: "110409A.wif"

Peak Name: "PFOA" Mass(es): "413.0/369.0 amu"


Comment: "Annotation: "



Sample Name: "SS33617" Sample ID: "Calibration Standard, 0.2 ng/mL" File: "110409A.wif"

Peak Name: "13CPFOA(M+2)(IS)" Mass(es): "415.0/370.0 amu"


Comment: "Annotation: "



Sample Name: "SS33616" Sample ID: "Calibration Standard, 0.5 ng/mL" File: "110409A.wif"

Peak Name: "PFOA" Mass(es): "413.0/369.0 amu"


Comment: "Annotation: "



Sample Name: "SS33616" Sample ID: "Calibration Standard, 0.5 ng/mL" File: "110409A.wif"

Peak Name: "13CPFOA(M+2)(IS)" Mass(es): "415.0/370.0 amu"

Comment: "Annotation: "









I Research, Inc.

Initials: MND

alyst Version: 1.4.2
rinting Time: 8:30:11 AM
rinting Date: Thursday, November 05, 2009

MPI Study: L19346
MPI Set No.: 110409A

Operator: Mark Neeley
Instrument No.: LC/MS/MS #9



alyst Version: 1.4.2
rinting Time: 8:30:12 AM

MPI Study: L19346
MPI Set No.: 110409A

Operator: Mark Neeley
Instrument No.: LC/MS/MS #9

rinting Date: Thursday, November 05, 2009

Sample Name: "Methanol Wash" Sample ID: "Methanol Wash" File: "110409A.wif"

Peak Name: "PFOA" Mass(es): "413.0/369.0 amu"

Conc. = 0.00 ng/mL
Calculated Conc: N/A
Acq. Date: 11/4/2009
Acq. Time: 5:00:06 PM

Modified: No

Proc. Algorithm: Analyst Classic

Bunching Factor: 1

Noise Threshold: 40.01 cps

Area Threshold: 200.04 cps

Smooths: 1

Width: 0.20

Height: 0.01

Peak Ratio: 5.00

Exp. Adj. Ratio: 4.00

Exp. Val. Ratio: 3.00

RT Window: 30.0 sec

Expected RT: 3.99 min

Use Relative RT: No


Int. Type: Base To Base

Retention Time: 3.94 min

Height: 1476 counts

Start Time: 3.83 min

End Time: 4.00 min



Sample Name: "Methanol Wash" Sample ID: "Methanol Wash" File: "110409A.wif"

Peak Name: "13CPFOA(M+2)(S)" Mass(es): "415.0/370.0 amu"

Conc. = 0.00 ng/mL
Calculated Conc: N/A
Acq. Date: 11/4/2009
Acq. Time: 5:00:06 PM

Modified: No

Proc. Algorithm: Analyst Classic

Bunching Factor: 1

Noise Threshold: 40.01 cps

Area Threshold: 200.04 cps

Smooths: 1

Width: 0.20

Height: 0.01

Peak Ratio: 5.00

Exp. Adj. Ratio: 4.00

Exp. Val. Ratio: 3.00

RT Window: 30.0 sec

Expected RT: 3.99 min

Use Relative RT: No


Int. Type: Base To Base

Retention Time: 3.94 min

Height: 1476 counts

Start Time: 3.83 min

End Time: 4.00 min



Sample Name: "Methanol Wash" Sample ID: "Methanol Wash" File: "110409A.wif"

Peak Name: "PFOA" Mass(es): "413.0/369.0 amu"

Conc. = 0.00 ng/mL
Calculated Conc: N/A
Acq. Date: 11/4/2009
Acq. Time: 5:11:01 PM

Modified: No

Proc. Algorithm: Analyst Classic

Bunching Factor: 1

Noise Threshold: 40.01 cps

Area Threshold: 200.04 cps

Smooths: 1

Width: 0.20

Sep. Height: 0.01

Exp. Peak Ratio: 5.00

Exp. Adj. Ratio: 4.00

Exp. Val. Ratio: 3.00

RT Window: 30.0 sec

Expected RT: 3.99 min

Use Relative RT: No


Int. Type: Base To Base

Retention Time: 3.93 min

Height: 1517 counts

Start Time: 3.82 min

End Time: 4.11 min



Sample Name: "Methanol Wash" Sample ID: "Methanol Wash" File: "110409A.wif"

Peak Name: "13CPFOA(M+2)(S)" Mass(es): "415.0/370.0 amu"

Conc. = 0.00 ng/mL
Calculated Conc: N/A
Acq. Date: 11/4/2009
Acq. Time: 5:11:01 PM

Modified: No

Proc. Algorithm: Analyst Classic

Bunching Factor: 1

Noise Threshold: 40.01 cps

Area Threshold: 200.04 cps

Smooths: 1

Width: 0.20

Sep. Height: 0.01

Exp. Peak Ratio: 5.00

Exp. Adj. Ratio: 4.00

Exp. Val. Ratio: 3.00

RT Window: 30.0 sec

Expected RT: 3.99 min

Use Relative RT: No


Int. Type: Base To Base

Retention Time: 3.94 min

Height: 1476 counts

Start Time: 3.83 min

End Time: 4.00 min



Sample Name: "Control" Sample ID: "MC4314 Deer Muscle Control" File: "110409A.wif"

Peak Name: "PFOA" Mass(es): "413.0/369.0 amu"

Conc. = 0.00 ng/mL
Calculated Conc: N/A
Acq. Date: 11/4/2009
Acq. Time: 5:21:55 PM

Modified: No

Proc. Algorithm: Analyst Classic

Bunching Factor: 1

Noise Threshold: 40.01 cps

Area Threshold: 200.04 cps

Smooths: 1

Width: 0.20

Sep. Height: 0.01

Exp. Peak Ratio: 5.00

Exp. Adj. Ratio: 4.00

Exp. Val. Ratio: 3.00

RT Window: 30.0 sec

Expected RT: 3.99 min

Use Relative RT: No


Int. Type: Base To Base

Retention Time: 3.93 min

Height: 12498 counts

Start Time: 3.82 min

End Time: 4.02 min



Sample Name: "Control" Sample ID: "MC4314 Deer Muscle Control" File: "110409A.wif"

Peak Name: "13CPFOA(M+2)(S)" Mass(es): "415.0/370.0 amu"

Conc. = 1.00 ng/mL
Calculated Conc: N/A
Acq. Date: 11/4/2009
Acq. Time: 5:21:55 PM

Modified: No

Proc. Algorithm: Analyst Classic

Bunching Factor: 1

Noise Threshold: 256.47 cps

Area Threshold: 1282.36 cps

Smooths: 1

Width: 0.20

Sep. Height: 0.01

Exp. Peak Ratio: 5.00

Exp. Adj. Ratio: 4.00

Exp. Val. Ratio: 1.00

RT Window: 30.0 sec

Expected RT: 3.99 min

Use Relative RT: No


Int. Type: Base To Base

Retention Time: 3.93 min

Height: 630414 counts

Start Time: 3.82 min

End Time: 4.13 min



alyst Version: 1.4.2
 rinting Time: 8:30:12 AM
 rinting Date: Thursday, November 05, 2009

MPI Study: L19346
 MPI Set No.: 110409A

Operator: Mark Neeley
 Instrument No.: LC/MS/MS #9

Sample Name: "LCS A" Sample ID: "MC4314 Deer Muscle Spike A, 10 ng/mL" File: "110409A.wif"

Peak Name: "PFOA" Mass(es): "413.0/369.0 amu"

Comment: " Annotation: "

Sample Index: 10

Sample Type: QC

Concentration: 10.0 ng/mL

Calculated Conc: 9.36 ng/mL

Acq. Date: 11/4/2009

Acq. Time: 5:32:53 PM

RT Window: 10.0 sec

Expected RT: 3.99 min

Relative RT: No

Int. Type: Base To Base

Retention Time: 3.93 min

Area: 311604 counts

Height: 114798 cps

Start Time: 3.87 min

End Time: 4.11 min

Intensity, cps

5.00e4

2.50e4

2.00e4

1.50e4

1.00e4

5.00e3

0.00

Time, min

Sample Name: "LCS A" Sample ID: "MC4314 Deer Muscle Spike A, 10 ng/mL" File: "110409A.wif"

Peak Name: "13CPFOA(M+2)(S)" Mass(es): "415.0/370.0 amu"

Comment: " Annotation: "

Sample Index: 10

Sample Type: QC

Concentration: 1.00 ng/mL

Calculated Conc: N/A

Acq. Date: 11/4/2009

Acq. Time: 5:32:53 PM

RT Window: 10.0 sec

Expected RT: 3.98 min

Relative RT: No

Int. Type: Base To Base

Retention Time: 3.93 min

Area: 765110 counts

Height: 284507 cps

Start Time: 3.87 min

End Time: 4.14 min

Intensity, cps

2.0e5

1.8e5

1.6e5

1.4e5

1.2e5

1.0e5

8.0e4

6.0e4

4.0e4

2.0e4

0.0

Time, min

Sample Name: "LCS B" Sample ID: "MC4314 Deer Muscle Spike B, 50 ng/mL" File: "110409A.wif"

Peak Name: "PFOA" Mass(es): "413.0/369.0 amu"

Comment: " Annotation: "

Sample Index: 11

Sample Type: QC

Concentration: 50.0 ng/mL

Calculated Conc: 47.1 ng/mL

Acq. Date: 11/4/2009

Acq. Time: 5:43:51 PM

RT Window: 30.0 sec

Expected RT: 3.99 min

Relative RT: No

Int. Type: Base To Base

Retention Time: 3.93 min

Area: 1418119 counts

Height: 520421 cps

Start Time: 3.87 min

End Time: 4.22 min

Intensity, cps

2.5e5

1.5e5

1.0e5

5.0e4

0.0

Time, min

Sample Name: "LCS B" Sample ID: "MC4314 Deer Muscle Spike B, 50 ng/mL" File: "110409A.wif"

Peak Name: "13CPFOA(M+2)(S)" Mass(es): "415.0/370.0 amu"

Comment: " Annotation: "

Sample Index: 11

Sample Type: QC

Concentration: 1.00 ng/mL

Calculated Conc: N/A

Acq. Date: 11/4/2009

Acq. Time: 5:43:51 PM

RT Window: 30.0 sec

Expected RT: 3.98 min

Relative RT: No

Int. Type: Base To Base

Retention Time: 3.92 min

Area: 709162 counts

Height: 263777 cps

Start Time: 3.87 min

End Time: 4.13 min

Intensity, cps

1.2e5

1.0e5

8.0e4

6.0e4

4.0e4

2.0e4

0.0

Time, min

Sample Name: "L19346-3 Spk C" Sample ID: "Deer #7 3.5 yr male-muscle Spike C, 50 ng/mL" File: "110409A.wif"

Peak Name: "PFOA" Mass(es): "413.0/369.0 amu"

Comment: " Annotation: "

Sample Index: 12

Sample Type: QC

Concentration: 50.0 ng/mL

Calculated Conc: 48.5 ng/mL

Acq. Date: 11/4/2009

Acq. Time: 5:54:48 PM

RT Window: 30.0 sec

Expected RT: 3.99 min

Relative RT: No

Int. Type: Base To Base

Retention Time: 3.93 min

Area: 1437886 counts

Height: 525845 cps

Start Time: 3.87 min

End Time: 4.18 min

Intensity, cps

2.5e5

1.5e5

1.0e5

5.0e4

0.0

Time, min

Sample Name: "L19346-3 Spk C" Sample ID: "Deer #7 3.5 yr male-muscle Spike C, 50 ng/mL" File: "110409A.wif"

Peak Name: "13CPFOA(M+2)(S)" Mass(es): "415.0/370.0 amu"

Comment: " Annotation: "

Sample Index: 12

Sample Type: QC

Concentration: 1.00 ng/mL

Calculated Conc: N/A

Acq. Date: 11/4/2009

Acq. Time: 5:54:48 PM

RT Window: 30.0 sec

Expected RT: 3.98 min

Relative RT: No

Int. Type: Base To Base

Retention Time: 3.93 min

Area: 699442 counts

Height: 260425 cps

Start Time: 3.87 min

End Time: 4.13 min

Intensity, cps

1.2e5

1.0e5

8.0e4

6.0e4

4.0e4

2.0e4

0.0

Time, min

alyst Version: 1.4.2

rinting Time: 8:30:12 AM

MPI Study: L19346

MPI Set No.: 110409A

Operator: Mark Neeley

Instrument No.: LC/MS/MS #9

rinting Date: Thursday, November 05, 2009

Sample Name: "L19346-1 Spk D" Sample ID: "Deer # 6 0.5 yr female-muscle Spike D, 50 ng/mL" File: "110409A.wif"


Peak Name: "PFDA" Mass(es): "413.0/369.0 amu"

Comment: "Annotation: "

Sample Index: 13
 Sample Type: QC
 Concentration: 50.0 ng/mL
 Calculated Conc: 46.3 ng/mL
 Acq. Date: 11/4/2009
 Acq. Time: 6:05:45 PM

Modified: No
 Proc. Algorithm: Analyst Classic
 Bunching Factor: 1
 Noise Threshold: 40.01 cps
 Area Threshold: 200.04 cps
 1. Smooths: 1
 2. Width: 0.20
 3. Height: 0.01
 4. Peak Ratio: 5.00
 5. Adj. Ratio: 4.00
 6. Exp. Val. Ratio: 3.00
 RT Window: 30.0 sec
 Expected RT: 3.99 min
 Use Relative RT: No

Intensity(cps)



Sample Name: "L19346-1 Spk D" Sample ID: "Deer # 6 0.5 yr female-muscle Spike D, 50 ng/mL" File: "110409A.wif"

Peak Name: "13CPFOA(M+2)(S)" Mass(es): "415.0/370.0 amu"

Comment: "Annotation: "

Sample Index: 13
 Sample Type: QC
 Concentration: 1.00 ng/mL
 Calculated Conc: N/A
 Acq. Date: 11/4/2009
 Acq. Time: 6:05:45 PM

Modified: No

Proc. Algorithm: Analyst Classic

Bunching Factor: 1

Noise Threshold: 256.47 cps

Area Threshold: 1282.36 cps

Num. Smooths: 1

Sep. Width: 0.20

Sep. Height: 0.01

Exp. Peak Ratio: 5.00

Exp. Adj. Ratio: 4.00

Exp. Val. Ratio: 3.00

RT Window: 30.0 sec

Expected RT: 3.98 min

Use Relative RT: No

Intensity(cps)



Sample Name: "L19346-3" Sample ID: "Deer # 7 3.5 yr male-muscle" File: "110409A.wif"


Peak Name: "PFDA" Mass(es): "413.0/369.0 amu"

Comment: "Annotation: "

Sample Index: 14
 Sample Type: Unknown
 Concentration: < 0
 Calculated Conc: N/A
 Acq. Date: 11/4/2009
 Acq. Time: 6:16:42 PM

Modified: No
 Proc. Algorithm: Analyst Classic
 Bunching Factor: 1
 Noise Threshold: 40.01 cps
 Area Threshold: 200.04 cps
 1. Smooths:
 2. Width: 0.20
 3. Height: 0.01
 4. Peak Ratio: 5.00
 5. Adj. Ratio: 4.00
 6. Exp. Val. Ratio: 3.00
 RT Window: 30.0 sec
 Expected RT: 3.99 min
 Use Relative RT: No

Intensity(cps)



Sample Name: "L19346-3" Sample ID: "Deer # 7 3.5 yr male-muscle" File: "110409A.wif"


Peak Name: "13CPFOA(M+2)(S)" Mass(es): "415.0/370.0 amu"

Comment: "Annotation: "

Sample Index: 14
 Sample Type: Unknown
 Concentration: 1.00 ng/mL
 Calculated Conc: N/A
 Acq. Date: 11/4/2009
 Acq. Time: 6:16:42 PM

Modified: No
 Proc. Algorithm: Analyst Classic
 Bunching Factor: 1
 Noise Threshold: 256.47 cps
 Area Threshold: 1282.36 cps
 Num. Smooths: 1
 Sep. Width: 0.20
 Sep. Height: 0.01
 Exp. Peak Ratio: 5.00
 Exp. Adj. Ratio: 4.00
 Exp. Val. Ratio: 3.00
 RT Window: 30.0 sec
 Expected RT: 3.98 min
 Use Relative RT: No

Intensity(cps)



Sample Name: "L19346-3 Dup" Sample ID: "Deer # 7 3.5 yr male-muscle Duplicate" File: "110409A.wif"


Peak Name: "PFDA" Mass(es): "413.0/369.0 amu"

Comment: "Annotation: "

Sample Index: 15
 Sample Type: Unknown
 Concentration: N/A
 Calculated Conc: < 0
 Acq. Date: 11/4/2009
 Acq. Time: 6:27:39 PM

Modified: No
 Proc. Algorithm: Analyst Classic
 Bunching Factor: 1
 Noise Threshold: 40.01 cps
 Area Threshold: 200.04 cps
 Num. Smooths:
 Sep. Width: 0.20
 Sep. Height: 0.01
 Exp. Peak Ratio: 5.00
 Exp. Adj. Ratio: 4.00
 Exp. Val. Ratio: 3.00
 Window: 30.0 sec
 Expected RT: 3.99 min
 Use Relative RT: No

Intensity(cps)



Sample Name: "L19346-3 Dup" Sample ID: "Deer # 7 3.5 yr male-muscle Duplicate" File: "110409A.wif"

Peak Name: "13CPFOA(M+2)(S)" Mass(es): "415.0/370.0 amu"

Comment: "Annotation: "

Sample Index: 15
 Sample Type: Unknown
 Concentration: 1.00 ng/mL
 Calculated Conc: N/A
 Acq. Date: 11/4/2009
 Acq. Time: 6:27:39 PM

Modified: No
 Proc. Algorithm: Analyst Classic
 Bunching Factor: 1
 Noise Threshold: 256.47 cps
 Area Threshold: 1282.36 cps
 Num. Smooths: 1
 Sep. Width: 0.20
 Sep. Height: 0.01
 Exp. Peak Ratio: 5.00
 Exp. Adj. Ratio: 4.00
 Exp. Val. Ratio: 3.00
 RT Window: 30.0 sec
 Expected RT: 3.98 min
 Use Relative RT: No

Intensity(cps)



lyst Version: 1.4.2
rinting Time: 8:30:12 AM
rinting Date: Thursday, November 05, 2009

MPI Study: L19346
MPI Set No.: 110409A

Operator: Mark Neeley
Instrument No.: LC/MS/MS #9

Sample Name: "L19346-1" Sample ID: "Deer # 8.0 yr female-muscle" File: "110409A.wif"

Peak Name: "PFOA" Mass(es): "413.0/369.0 amu"

Comment: " Annotation: "

Sample Index: 16

Sample Type: Unknown

Concentration: N/A

Calculated Conc: < 0

Acq. Date: 11/4/2009

Acq. Time: 6:38:37 PM

Modified: No

Proc. Algorithm: Analyst Classic

Bunching Factor: 1

Noise Threshold: 40.01 cps

Area Threshold: 200.04 cps

Smooths: 1

Sep. Width: 0.20

Sep. Height: 0.01

Exp. Peak Ratio: 5.00


Exp. Adj. Ratio: 4.00

Exp. Val. Ratio: 3.00

RT Window: 30.0 sec

Expected RT: 3.99 min

Relative RT: No



Sample Name: "L19346-1" Sample ID: "Deer # 8.0 yr female-muscle" File: "110409A.wif"

Peak Name: "13CPFOA(M+2)(IS)" Mass(es): "415.0/370.0 amu"

Comment: " Annotation: "

Sample Index: 16

Sample Type: Unknown

Concentration: 1.00 ng/mL

Calculated Conc: N/A

Acq. Date: 11/4/2009

Acq. Time: 6:38:37 PM

Modified: No

Proc. Algorithm: Analyst Classic

Bunching Factor: 1

Noise Threshold: 256.47 cps

Area Threshold: 1282.36 cps

Num. Smooths: 1

Sep. Width: 0.20

Sep. Height: 0.01

Exp. Peak Ratio: 5.00


Exp. Adj. Ratio: 4.00

Exp. Val. Ratio: 3.00

RT Window: 30.0 sec

Expected RT: 3.98 min

Relative RT: No



Sample Name: "L19346-1 Dup" Sample ID: "Deer # 8.0 yr female-muscle Duplicate" File: "110409A.wif"

Peak Name: "PFOA" Mass(es): "413.0/369.0 amu"

Comment: " Annotation: "

Sample Index: 17

Sample Type: Unknown

Concentration: N/A

Calculated Conc: < 0

Acq. Date: 11/4/2009

Acq. Time: 6:49:35 PM

Modified: Yes

Proc. Algorithm: Analyst Classic

Bunching Factor: 1

Noise Threshold: 40.01 cps

Area Threshold: 200.04 cps

Smooths: 1

Sep. Width: 0.20

Sep. Height: 0.01

Exp. Peak Ratio: 5.00


Exp. Adj. Ratio: 4.00

Exp. Val. Ratio: 3.00

RT Window: 30.0 sec

Expected RT: 3.92 min

Relative RT: No



Sample Name: "L19346-1 Dup" Sample ID: "Deer # 8.0 yr female-muscle Duplicate" File: "110409A.wif"

Peak Name: "13CPFOA(M+2)(IS)" Mass(es): "415.0/370.0 amu"

Comment: " Annotation: "

Sample Index: 17

Sample Type: Unknown

Concentration: 1.00 ng/mL

Calculated Conc: N/A

Acq. Date: 11/4/2009

Acq. Time: 6:49:35 PM

Modified: Yes

Proc. Algorithm: Analyst Classic

Bunching Factor: 1

Noise Threshold: 256.47 cps

Area Threshold: 1282.36 cps

Num. Smooths: 1

Sep. Width: 0.20

Sep. Height: 0.01

Exp. Peak Ratio: 5.00


Exp. Adj. Ratio: 4.00

Exp. Val. Ratio: 3.00

RT Window: 30.0 sec

Expected RT: 3.98 min

Relative RT: No



Sample Name: "SS33815" Sample ID: "CCV, 1.0 ng/mL" File: "110409A.wif"

Peak Name: "PFOA" Mass(es): "413.0/369.0 amu"

Comment: " Annotation: "

Sample Index: 18

Sample Type: QC

Concentration: 1.00 ng/mL

Calculated Conc: 1.04 ng/mL

Acq. Date: 11/4/2009

Acq. Time: 7:00:33 PM

Modified: No

Proc. Algorithm: Analyst Classic

Bunching Factor: 1

Noise Threshold: 40.01 cps

Area Threshold: 200.04 cps

Smooths: 1

Sep. Width: 0.20

Sep. Height: 0.03

Exp. Peak Ratio: 5.00


Exp. Adj. Ratio: 4.00

Exp. Val. Ratio: 3.00

RT Window: 30.0 sec

Expected RT: 3.99 min

Relative RT: No



Sample Name: "SS33815" Sample ID: "CCV, 1.0 ng/mL" File: "110409A.wif"

Peak Name: "13CPFOA(M+2)(IS)" Mass(es): "415.0/370.0 amu"

Comment: " Annotation: "

Sample Index: 18

Sample Type: QC

Concentration: 1.00 ng/mL

Calculated Conc: N/A

Acq. Date: 11/4/2009

Acq. Time: 7:00:33 PM

Modified: No

Proc. Algorithm: Analyst Classic

Bunching Factor: 1

Noise Threshold: 256.47 cps

Area Threshold: 1282.36 cps

Num. Smooths: 1

Sep. Width: 0.20

Sep. Height: 0.01

Exp. Peak Ratio: 5.00


Exp. Adj. Ratio: 4.00

Exp. Val. Ratio: 3.00

RT Window: 30.0 sec

Expected RT: 3.98 min

Relative RT: No




lyst Version: 1.4.2
rinting Time: 8:29:38 AM
rinting Date: Thursday, November 05, 2009

MPI Study: L19346
MPI Set No.: 110409A

Operator: Mark Neeley
Instrument No.: LC/MS/MS #9

Nov 11-5-09

09_110409A.rdb (PFOS): "Linear" Regression ("1 / x" weighting): $y = 1.07 x + 0.00957$ ($r = 0.9998$)









lyst Version: 1.4.2
rinting Time: 8:30:30 AM

MPI Study: L19346
MPI Set No.: 110409A

Operator: Mark Neeley
Instrument No.: LC/MS/MS #9

rinting Date: Thursday, November 05, 2009



Research, Inc.

Initials Mark

Page 1 of 6

Date 11-5-09

Sample

Index 1 To 16

lyst Version: 1.4.2
rinting Time: 8:30:30 AM
rinting Date: Thursday, November 05, 2009

MPI Study: L19346
MPI Set No.: 110409A

Operator: Mark Neeley
Instrument No.: LC/MS/MS #9

Sample Name: "SS33815" Sample ID: "Calibration Standard, 1.0 ng/mL" File: "110409A.wif"

Peak Name: "PFOS" Mass(es): "499.0/80.0 amu"

Comment: "Annotation: "

Sample Index: 4
Sample Type: Standard
Concentration: 1.00 ng/mL
Calculated Conc: 1.00 ng/mL
Acq. Date: 11/4/2009
Acq. Time: 4:27:18 PM
Modified: No
Proc. Algorithm: Analyst Classic
Bunching Factor: 1
Noise Threshold: 15.73 cps
Area Threshold: 79.65 cps
Smooths: 1
Sep. Width: 0.20
Sep. Height: 0.01
Exp. Peak Ratio: 5.00
Exp. Adj. Ratio: 4.00
Exp. Val. Ratio: 3.00
RT Window: 30.0 sec
Expected RT: 4.06 min
Relative RT: No
Type: Base To Base
Retention Time: 4.01 min
Area: 369109 counts
Height: 113187 cps
Start Time: 3.90 min
End Time: 4.25 min
3.00e4
2.50e4
2.00e4
1.50e4
1.00e4
5000.00
0.00
Intensity, cps
Time, min

Sample Name: "SS33814" Sample ID: "Calibration Standard, 2.0 ng/mL" File: "110409A.wif"

Peak Name: "PFOS" Mass(es): "499.0/80.0 amu"

Comment: "Annotation: "

Sample Index: 5
Sample Type: Standard
Concentration: 2.00 ng/mL
Calculated Conc: 2.07 ng/mL
Acq. Date: 11/4/2009
Acq. Time: 4:38:14 PM
Modified: No
Proc. Algorithm: Analyst Classic
Bunching Factor: 1
Noise Threshold: 15.73 cps
Area Threshold: 79.65 cps
Smooths: 1
Sep. Width: 0.20
Sep. Height: 0.01
Exp. Peak Ratio: 5.00
Exp. Adj. Ratio: 4.00
Exp. Val. Ratio: 3.00
RT Window: 30.0 sec
Expected RT: 4.06 min
Relative RT: No
Type: Base To Base
Retention Time: 4.01 min
Area: 747914 counts
Height: 235561 cps
Start Time: 3.91 min
End Time: 4.29 min
6.0e4
4.0e4
2.0e4
0.0
Intensity, cps
Time, min

Sample Name: "SS33813" Sample ID: "Calibration Standard, 5.0 ng/mL" File: "110409A.wif"

Peak Name: "PFOS" Mass(es): "499.0/80.0 amu"

Comment: "Annotation: "

Sample Index: 6
Sample Type: Standard
Concentration: 5.00 ng/mL
Calculated Conc: 4.93 ng/mL
Acq. Date: 11/4/2009
Acq. Time: 4:49:09 PM
Modified: No
Proc. Algorithm: Analyst Classic
Bunching Factor: 1
Noise Threshold: 15.73 cps
Area Threshold: 79.65 cps
Smooths: 1
Sep. Width: 0.20
Sep. Height: 0.01
Exp. Peak Ratio: 5.00
Exp. Adj. Ratio: 4.00
Exp. Val. Ratio: 3.00
RT Window: 30.0 sec
Expected RT: 4.06 min
Relative RT: No
Type: Base To Base
Retention Time: 4.00 min
Area: 1751107 counts
Height: 539144 cps
Start Time: 3.90 min
End Time: 4.29 min
2.5e5
1.0e5
5.0e4
0.0
Intensity, cps
Time, min

Sample Name: "SS33815" Sample ID: "Calibration Standard, 1.0 ng/mL" File: "110409A.wif"

Peak Name: "13CPFOS(M+4)IS" Mass(es): "503.0/80.0 amu"
Comment: "Annotation: "

Sample Index: 4
Sample Type: Standard
Concentration: 1.00 ng/mL
Calculated Conc: N/A 1.20e5
Acq. Date: 11/4/2009
Acq. Time: 4:27:18 PM
Modified: No
Proc. Algorithm: Analyst Classic
Bunching Factor: 1
Noise Threshold: 140.21 cps
Area Threshold: 701.07 cps
Smooths: 1
Sep. Width: 0.20
Sep. Height: 0.01
Exp. Peak Ratio: 5.00
Exp. Adj. Ratio: 4.00
Exp. Val. Ratio: 3.00
RT Window: 30.0 sec
Expected RT: 4.06 min
Relative RT: No
Type: Base To Base
Retention Time: 4.00 min
Area: 342687 counts
Height: 124540 cps
Start Time: 3.94 min
End Time: 4.20 min
3.00e4
2.00e4
1.00e4
0.00
Intensity, cps
Time, min

Sample Name: "SS33814" Sample ID: "Calibration Standard, 2.0 ng/mL" File: "110409A.wif"

Peak Name: "13CPFOS(M+4)IS" Mass(es): "503.0/80.0 amu"

Comment: "Annotation: "

Sample Index: 5
Sample Type: Standard
Concentration: 1.00 ng/mL
Calculated Conc: N/A 1.20e5
Acq. Date: 11/4/2009
Acq. Time: 4:38:14 PM
Modified: No
Proc. Algorithm: Analyst Classic
Bunching Factor: 1
Noise Threshold: 140.21 cps
Area Threshold: 701.07 cps
Smooths: 1
Sep. Width: 0.20
Sep. Height: 0.01
Exp. Peak Ratio: 5.00
Exp. Adj. Ratio: 4.00
Exp. Val. Ratio: 3.00
RT Window: 30.0 sec
Expected RT: 4.06 min
Relative RT: No
Type: Base To Base
Retention Time: 4.01 min
Area: 338050 counts
Height: 123612 cps
Start Time: 3.95 min
End Time: 4.21 min
3.00e4
2.00e4
1.00e4
0.00
Intensity, cps
Time, min

Sample Name: "SS33813" Sample ID: "Calibration Standard, 5.0 ng/mL" File: "110409A.wif"

Peak Name: "13CPFOS(M+4)IS" Mass(es): "503.0/80.0 amu"







Comment: "Annotation: "

Sample Index: 6
Sample Type: Standard
Concentration: 1.00 ng/mL
Calculated Conc: N/A 1.10e5
Acq. Date: 11/4/2009
Acq. Time: 4:49:09 PM
Modified: No
Proc. Algorithm: Analyst Classic
Bunching Factor: 1
Noise Threshold: 140.21 cps
Area Threshold: 701.07 cps
Smooths: 1
Sep. Width: 0.20
Sep. Height: 0.01
Exp. Peak Ratio: 5.00
Exp. Adj. Ratio: 4.00
Exp. Val. Ratio: 3.00
RT Window: 30.0 sec
Expected RT: 4.06 min
Relative RT: No
Type: Base To Base
Retention Time: 4.00 min
Area: 332437 counts
Height: 119161 cps
Start Time: 3.94 min
End Time: 4.21 min
3.00e4
2.00e4
1.00e4
0.00
Intensity, cps
Time, min

alyst Version: 1.4.2
rinting Time: 8:30:30 AM
rinting Date: Thursday, November 05, 2009

MPI Study: L19346
MPI Set No.: 110409A

Operator: Mark Neeley
Instrument No.: LC/MS/MS #9



ilyst Version: 1.4.2
Printing Time: 8:30:30 AM

MPI Study: L19346
MPI Set No.: 110409A


Operator: Mark Neeley
Instrument No.: LC/MS/MS #9

Printing Date: Thursday, November 05, 2009

Sample Name: "LCS A" Sample ID: "MC4314 Deer Muscle Spike A, 10 ng/mL" File: "110409A.wif"

Peak Name: "PFOS" Mass(es): "499.0/80.0 amu"


Comment: "Annotation: "



Sample Name: "LCS A" Sample ID: "MC4314 Deer Muscle Spike A, 10 ng/mL" File: "110409A.wif"

Peak Name: "13CPFOS(M+4)(IS)" Mass(es): "503.0/80.0 amu"


Comment: "Annotation: "



Sample Name: "LCS B" Sample ID: "MC4314 Deer Muscle Spike B, 50 ng/mL" File: "110409A.wif"

Peak Name: "PFOS" Mass(es): "499.0/80.0 amu"


Comment: "Annotation: "



Sample Name: "LCS B" Sample ID: "MC4314 Deer Muscle Spike B, 50 ng/mL" File: "110409A.wif"

Peak Name: "13CPFOS(M+4)(IS)" Mass(es): "503.0/80.0 amu"


Comment: "Annotation: "



Sample Name: "L19346-3 Spk C" Sample ID: "Deer #73.5 yr male-muscle Spike C, 50 ng/mL" File: "110409A.wif"

Peak Name: "PFOS" Mass(es): "499.0/80.0 amu"


Comment: "Annotation: "



Sample Name: "L19346-3 Spk C" Sample ID: "Deer #73.5 yr male-muscle Spike C, 50 ng/mL" File: "110409A.wif"

Peak Name: "13CPFOS(M+4)(IS)" Mass(es): "503.0/80.0 amu"

Comment: "Annotation: "









alyst Version: 1.4.2
rinting Time: 8:30:30 AM

MPI Study: L19346
MPI Set No.: 110409A

Operator: Mark Neeley
Instrument No.: LC/MS/MS #9







rinting Date: Thursday, November 05, 2009



alyst Version: 1.4.2
rinting Time: 8:30:30 AM
rinting Date: Thursday, November 05, 2009

MPI Study: L19346
MPI Set No.: 110409A

Operator: Mark Neeley
Instrument No.: LC/MS/MS #9



RAW DATA REPORT

Sponsor Study No: NA Limit of Quantitation: 10 ng/g Set No: 110409AD
 MPI Study No: L19345 Injection Volume: 15 µL Analyst: Mark Neely
 Analyte: PFOS Matrix: Deer Muscle Instrument Type: LC/MS/MS Unit # 9
 Ions Monitored: 499 -> 80 Extraction Date: 11/04/09
 Site: NA Analyzed on: 11/05/09

MPI Research ID	Sponsor ID	Sample Code	Sample Index No.	Internal				Internal		Amount		
				Std. Conc. (ng/mL)	Std. Conc. (ng/mL)	Aliquot Factor (AF)	Dilution Factor (DF)	Peak Area	Peak Area	Found (ng/g)	Analyte Added (ng/g)	Recovery (%)
SS33618	-	CS	1	0.100	1.0	-	-	37091	326346	-	-	-
SS33617	-	CS	2	0.200	1.0	-	-	71080	310927	-	-	-
SS33616	-	CS	3	0.500	1.0	-	-	172764	318796	-	-	-
SS33615	-	CS	4	1.00	1.0	-	-	344831	319636	-	-	-
SS33614	-	CS	5	2.00	1.0	-	-	702604	319828	-	-	-
SS33613	-	CS	6	5.00	1.0	-	-	1728169	318164	-	-	-
Methanol Wash	-	W	7	-	-	-	-	315	0	-	-	-
Methanol Wash	-	W	8	-	-	-	-	198	0	-	-	-
L19346-3 Spk C	Deer # 7 3.5 yr male-muscle Spike C	S	9	-	1.0	20	10	243599	345013	128	50	79
SS33615		CCV	10	1.00	1.0	-	-	365129	333224	1.00	1.0	100

Standard Curve: Linear (1/x weighted)

Intercept = 0.00598

Slope = 1.09

Coef. Of Det. = 1.00000

$$\text{Analyte Found (ng/g)} = (((\text{analyte peak area}/\text{IS peak area}) - \text{intercept}) / \text{slope}) \times \text{IS conc.} \times \text{AF} \times \text{DF}/\text{Sample weight}$$

$$\text{Recovery (\%)} = \frac{[\text{Analyte found (ng/g)} - \text{Analyte found in control (ng/g)}]}{\text{amount Analyte added (ng/g)}} \times 100$$

CS = Calibration standard

CCV = Continuing Calibration Verification

C = Control sample

S = Sample

LF = Lab fortified sample

FF = Field fortified sample

LCS = Laboratory Control Spike

W = Methanol Wash

ND = Not detected = Response between 0 and LOD

NQ = Not quantifiable = Response between LOD and LOQ

Spreadsheet prepared by: MON 111-6-09

alyst Version: 1.4.2
rinting Time: 2:00:54 PM

MPI Study: L19346
MPI Set No.: 110409AD

Operator: Mark Neeley
Instrument No.: LC/MS/MS #9

rinting Date: Thursday, November 05, 2009

MDW 11-5-09

object:\sc1wp5556\mdrive\PE SCIEX DATA\Projects\P5195 Batch:09_110409AD Muscle Tab:Sample Set:SET1 AcqMethod:P5195_10
Start

Sample Name	Sample ID	Vial Position	Data File
SS33618	Calibration Standard, 0.1 ng/mL	1	09_110409AD\110409AD
SS33617	Calibration Standard, 0.2 ng/mL	2	09_110409AD\110409AD
SS33616	Calibration Standard, 0.5 ng/mL	3	09_110409AD\110409AD
SS33615	Calibration Standard, 1.0 ng/mL	4	09_110409AD\110409AD
SS33614	Calibration Standard, 2.0 ng/mL	5	09_110409AD\110409AD
SS33613	Calibration Standard, 5.0 ng/mL	6	09_110409AD\110409AD
Methanol Wash	Methanol Wash	91	09_110409AD\110409AD
Methanol Wash	Methanol Wash	91	09_110409AD\110409AD
L19346-3 Spk C	Deer # 7 3.5 yr male-muscle Spike C, 50 ng/mL, DF=50	90	09_110409AD\110409AD
SS33615	CCV, 1.0 ng/mL	4 <i>11-5-09</i>	09_110409AD\110409AD

LC/MS/MS SYSTEM AND OPERATING CONDITIONS

Protocol No: NA

MPI Study No: L19346

Instrument: AB API 4000 Biomolecular Mass Analyzer, (LC/MS/MS #9)
SCIEX Turbo Ion Spray Liquid Introduction Interface
Turbo Ion spray temperature = 450 °C

Computer: Dell OptiPlex GX 110

Software: PE Sciex Analyst 1.4

HPLC Equipment: Hewlett Packard (HP) Series 1100
HP Quat Pump HP Vacuum Degasser
HP Autosampler HP Column Oven

HPLC Column: Phenomenex Luna C8 (2) Mercury, 2cm x 4mm, 3 µm (ExyLIMS ID:
MA0052622)

Column Temperature: 35°C

Mobile Phase (A): 2 mM Ammonium Acetate in Water (ExyLIMS ID: SL0045925)

Mobile Phase (B): Methanol (ExyLIMS ID: RE0047880)

Injected Volume: 15 µL

<u>Time (min)</u>	<u>% A</u>	<u>% B</u>	<u>Flow Rate (µL/min)</u>
0.0	90	10	750
0.5	90	10	750
2.0	10	90	750
5.0	10	90	750
5.1	0	100	750
6.0	0	100	750
6.1	90	10	750
10.0	90	10	750

Ions monitored:

<u>Analyte</u>	<u>Parent ion</u>	<u>Daughter ion(s)</u>	<u>Dwell (secs)</u>
PFOA	413	369	0.200
PFOS	499	80	0.200
¹³ C PFOA (m+2)	415	370	0.200
Internal Standard			
¹³ C PFOS (m+4)	503	80	0.200
Internal Standard			

Analyst:

Mark Neeley *11-5-09*

MPI Research, Inc.

3058 Research Drive, State College, PA 16801

Phone: (814) 272-1039 FAX: (814) 231-1580

All Handwritten Peak ID's by: Mark 11-6-09

Analyst Version: 1.4.2
Printing Time: 2:01:15 PM
Printing Date: Thursday, November 05, 2009

MPI Study: L19346
MPI Set No.: 110409AD

Operator: Mark Neeley
Instrument No.: LC/MS/MS #9

MDH 11-5-09

Acquisition Information:

Acquisition Method: P5195_102909.dam
Created: Thursday October 29 2009 09: 35: 12 AM
Last Modified: Friday October 30 2009 10: 56: 20 AM
Comment: PFQA PFOS
Synchronization Mode: LC Sync
Auto-Equilibration: Off
Acquisition Duration: 10min0sec
Number Of Scans: 732
Periods In File: 1
Acquisition Module: Acquisition Method
Software version: Analyst 1.4.2

Period 1:

Scans in Period: 732
Relative Start Time: 0.00 msec
Experiments in Period: 1

Period 1 Experiment 1:

Scan Type: MRM (MRM)
Polarity: Negative
Scan Mode: N/A
Ion Source: Turbo Spray
Resolution Q1: Unit
Resolution Q3: Unit
Intensity Thres.: 0.00 cps
Settling Time: 0.0000 msec
MR Pause: 5.0070 msec
MCA: No
Step Size: 0.00 amu

Q1 Mass (amu)	Q3 Mass (amu)	Dwell (msec)	Param	Start	Stop
413.00	369.00	200.00	DF	-32.00	-32.00
			CE	-18.00	-18.00
			CXE	-13.20	-13.20

Q1 Mass (amu)	Q3 Mass (amu)	Dwell (msec)	Param	Start	Stop
415.00	370.00	200.00	DF	-32.00	-32.00
			CE	-18.00	-18.00
			CXE	-13.20	-13.20

Q1 Mass (amu)	Q3 Mass (amu)	Dwell (msec)	Param	Start	Stop
499.00	80.00	200.00	DF	-83.00	-83.00
			CE	-88.00	-88.00
			CXE	-6.00	-6.00

Q1 Mass (amu)	Q3 Mass (amu)	Dwell (msec)	Param	Start	Stop
500.00	80.00	200.00	DF	-83.00	-83.00
			CE	-88.00	-88.00

Clyst Version: 1.4.2
rinting Time: 2:01:15 PM
rinting Date: Thursday, November 05, 2009

MPI Study: L19346
MPI Set No.: 110409AD

Operator: Mark Neeley
Instrument No.: LC/MS/MS #9

MDN 11-5-09

CXF -6.00 -6.00

Parameter Table (Period 1 Experiment 1):

AD: 7.00
CUR: 20.00
S1: 50.00
S2: 40.00
S: -4500.00
TEM: 450.00
He: ON
F: -10.00

Agilent LC Pump Method Properties

Pump Model: Agilent 1100 LC Quaternary Pump
Minimum Pressure (psi): 0.0
Maximum Pressure (psi): 5801.0
Compressibility: 100.0
Dead Volume (μ l): 40.0
Stroke Volume (μ l): -1.0
Maximum Flow Ramp (ml/min²): 100.0
Maximum Pressure Ramp (psi/sec): 290.0

Step Table:

Step	Total Time (min)	Flow Rate (μ l/min)	A (%)	B (%)	C (%)	D (%)	TE#1	TE#2	TE#3	TE#4
0	0.00	750	0.0	0.0	90.0	10.0	open	open	open	open
1	0.50	750	0.0	0.0	90.0	10.0	open	open	open	open
2	2.00	750	0.0	0.0	10.0	90.0	open	open	open	open
3	5.00	750	0.0	0.0	10.0	90.0	open	open	open	open
4	5.10	750	0.0	0.0	0.0	100.0	open	open	open	open
5	6.00	750	0.0	0.0	0.0	100.0	open	open	open	open
6	6.10	750	0.0	0.0	90.0	10.0	open	open	open	open
7	10.00	750	0.0	0.0	90.0	10.0	open	open	open	open

Primary Flow Rate (μ l/min): 200.0

Flow Sensor Calibration Table Index: 0

Agilent Column Oven Properties

Left Temperature (°C): 35.00

Right Temperature (°C): 35.00

Temperature Tolerance +/- (°C): 1.00

Start Acquisition Tolerance +/- (°C): 0.50

Time Table (Not Used)

Column Switching Valve Installed

Position for first sample in the batch: Left

Use same position for all samples in the batch

Agilent Autosampler Properties

Autosampler Model: Agilent 1100 Wellplate Autosampler

Syringe Size (μ l): 100

Injection Volume (μ l): 15.00

Draw Speed (μ l/min): 200.0

Eject Speed (μ l/min): 200.0

Needle Level (mm): 0.00

Temperature Control Enabled

Setpoint (4 - 40 C): 4

Wash Location: Flush Port

Wash Time (1 - 999 sec): 10

Automatic Delay Volume Reduction Not Used

PI Research, Inc.

lyst Version: 1.4.2
rinting Time: 2:01:16 PM
rinting Date: Thursday, November 05, 2009

MPI Study: L19346
MPI Set No.: 110409AD

Operator: Mark Neeley
Instrument No.: LC/MS/MS #9

MDW 11-5-09

Equilibration Time (sec): 2
Enable Vial/Well Bottom Sensing No
se Custom Injector Program No

Analyst Version: 1.4.2
Printing Time: 8:43:43 AM
Printing Date: Friday, November 06, 2009


MPI Study: L19346
MPI Set No.: 110409AD

Operator: Mark Neeley
Instrument No.: LC/MS/MS #9

MDP 11-6-09

09_110409AD.rdb (PFOS): "Linear" Regression ("1 / x" weighting): $y = 1.09 x + 0.00598$ ($r = 1.0000$)

Analyte Area / IS Area



lyst Version: 1.4.2
rinting Time: 8:44:29 AM
rinting Date: Friday, November 06, 2009

MPI Study: L19346
MPI Set No.: 110409AD

Operator: Mark Neeley
Instrument No.: LC/MS/MS #9

Sample Name: "SS33615" Sample ID: "Calibration Standard, 1.0 ng/mL" File: "110409AD.wif"

Peak Name: "PFOS" Mass(es): "499.0/80.0 amu"

Comment: " Annotation: "

Sample Index: 4

Sample Type: Standard

Concentration: 1.00 ng/mL

Calculated Conc: 0.988 ng/mL

Acq. Date: 11/5/2009

Acq. Time: 3:58:25 PM

Modified: No

Proc. Algorithm: Analyst Classic

Bunching Factor: 1

Noise Threshold: 15.73 cps

Area Threshold: 78.65 cps

Num. Smooths: 1

Sep. Width: 0.20

Sep. Height: 0.01

Exp. Peak Ratio: 5.00

Exp. Adj. Ratio: 4.00

Exp. Val. Ratio: 3.00

RT Window: 30.0 sec

Expected RT: 4.06 min

Relative RT: No

Int. Type: Base To Base

Retention Time: 4.01 min

Area: 344831 counts

Height: 102159 cps

Start Time: 3.91 min

End Time: 4.26 min

Intensity, cps

4.01

Time min

0.0 2 4 6 8

Sample Name: "SS33614" Sample ID: "Calibration Standard, 2.0 ng/mL" File: "110409AD.wif"

Peak Name: "PFOS" Mass(es): "499.0/80.0 amu"

Comment: " Annotation: "

Sample Index: 5

Sample Type: Standard

Concentration: 2.00 ng/mL

Calculated Conc: 2.02 ng/mL

Acq. Date: 11/5/2009

Acq. Time: 4:09:19 PM

Modified: No

Proc. Algorithm: Analyst Classic

Bunching Factor: 1

Noise Threshold: 15.73 cps

Area Threshold: 78.65 cps

Num. Smooths: 1

Sep. Width: 0.20

Sep. Height: 0.01

Exp. Peak Ratio: 5.00

Exp. Adj. Ratio: 4.00

Exp. Val. Ratio: 3.00

RT Window: 30.0 sec

Expected RT: 4.06 min

Relative RT: No

Int. Type: Base To Base

Retention Time: 4.01 min

Area: 702604 counts

Height: 210047 cps

Start Time: 3.91 min

End Time: 4.39 min

Intensity, cps

4.01

Time min

0.0 2 4 6 8

Sample Name: "SS33613" Sample ID: "Calibration Standard, 5.0 ng/mL" File: "110409AD.wif"

Peak Name: "PFOS" Mass(es): "499.0/80.0 amu"

Comment: " Annotation: "

Sample Index: 6

Sample Type: Standard

Concentration: 5.00 ng/mL

Calculated Conc: 5.00 ng/mL

Acq. Date: 11/5/2009

Acq. Time: 4:20:14 PM

Modified: No

Proc. Algorithm: Analyst Classic

Bunching Factor: 1

Noise Threshold: 15.73 cps

Area Threshold: 78.65 cps

Num. Smooths: 1

Sep. Width: 0.20

Sep. Height: 0.01

Exp. Peak Ratio: 5.00

Exp. Adj. Ratio: 4.00

Exp. Val. Ratio: 3.00

RT Window: 30.0 sec

Expected RT: 4.06 min

Relative RT: No

Int. Type: Base To Base

Retention Time: 4.01 min

Area: 1728169 counts

Height: 514006 cps

Start Time: 3.90 min

End Time: 4.55 min

Intensity, cps

4.01

Time min

0.0 2 4 6 8

PI Research, Inc.

Sample Name: "Methanol Wash" Sample ID: "Methanol Wash" File: "110409AD.wiff"

Peak Name: "PFOS" Mass(es): "499.0/60.0 amu"

Comment: "Annotation: "

Sample Index: 7

Sample Type: Solvent

Concentration: 0.00 ng/mL

Calculated Conc: N/A

Acq. Date: 11/5/2009

Acq. Time: 4:31:09 PM

Modified: No

Proc. Algorithm: Analyst Classic

Bunching Factor: 1

Noise Threshold: 15.73 cps

Area Threshold: 78.65 cps

Smooths: 1

Sep. Width: 0.20

Sep. Height: 0.01

Exp. Peak Ratio: 5.00

Exp. Adj. Ratio: 4.00

Exp. Val. Ratio: 3.00

RT Window: 30.0 sec

Expected RT: 4.06 min

Use Relative RT: No

Int. Type: Base To Base


Retention Time: 4.00 min

Area: 315 counts

Height: 71 cps

Start Time: 3.92 min

End Time: 4.06 min



Sample Name: "Methanol Wash" Sample ID: "Methanol Wash" File: "110409AD.wiff"

Peak Name: "13CPFOS(M+4)(S)" Mass(es): "503.0/60.0 amu"

Comment: "Annotation: "

Sample Index: 7

Sample Type: Solvent

Concentration: 1.00 ng/mL

Calculated Conc: N/A

Acq. Date: 11/5/2009

Acq. Time: 4:31:09 PM

Modified: No

Proc. Algorithm: Analyst Classic

Bunching Factor: 1

Noise Threshold: 15.73 cps

Area Threshold: 78.65 cps

Smooths: 1

Sep. Width: 0.20

Sep. Height: 0.01

Exp. Peak Ratio: 5.00

Exp. Adj. Ratio: 4.00

Exp. Val. Ratio: 3.00

RT Window: 30.0 sec

Expected RT: 4.06 min

Use Relative RT: No

Int. Type: Base To Base


Retention Time: 4.00 min

Area: 0.27 counts

Height: 225 cps

Start Time: 3.92 min

End Time: 4.06 min



Sample Name: "Methanol Wash" Sample ID: "Methanol Wash" File: "110409AD.wiff"

Peak Name: "PFOS" Mass(es): "499.0/60.0 amu"

Comment: "Annotation: "

Sample Index: 9

Sample Type: Solvent

Concentration: 0.00 ng/mL

Calculated Conc: N/A

Acq. Date: 11/5/2009

Acq. Time: 4:42:02 PM

Modified: No

Proc. Algorithm: Analyst Classic

Bunching Factor: 1

Noise Threshold: 15.73 cps

Area Threshold: 78.65 cps

Smooths: 1

Sep. Width: 0.20

Sep. Height: 0.01

Exp. Peak Ratio: 5.00

Exp. Adj. Ratio: 4.00

Exp. Val. Ratio: 3.00

RT Window: 30.0 sec

Expected RT: 4.06 min

Use Relative RT: No

Int. Type: Base To Base


Retention Time: 4.00 min

Area: 198 counts

Height: 71 cps

Start Time: 3.92 min

End Time: 4.07 min



Sample Name: "Methanol Wash" Sample ID: "Methanol Wash" File: "110409AD.wiff"

Peak Name: "13CPFOS(M+4)(S)" Mass(es): "503.0/60.0 amu"

Comment: "Annotation: "

Sample Index: 9

Sample Type: Solvent

Concentration: 1.00 ng/mL

Calculated Conc: N/A

Acq. Date: 11/5/2009

Acq. Time: 4:42:02 PM

Modified: No

Proc. Algorithm: Analyst Classic

Bunching Factor: 1

Noise Threshold: 15.73 cps

Area Threshold: 78.65 cps

Smooths: 1

Sep. Width: 0.20

Sep. Height: 0.01

Exp. Peak Ratio: 5.00

Exp. Adj. Ratio: 4.00

Exp. Val. Ratio: 3.00

RT Window: 30.0 sec

Expected RT: 4.06 min

Use Relative RT: No

Int. Type: Base To Base


Retention Time: 4.00 min

Area: 0.28 counts

Height: 3.29 cps

Start Time: 4.00 min

End Time: 5.48 min



Sample Name: "L19348-3 Spk C" Sample ID: "Deer #7 3.5 yr male-muscle Spike C, 50 ng/mL, DF=10" File: "110409AD.wiff"

Peak Name: "PFOS" Mass(es): "499.0/60.0 amu"

Comment: "Annotation: "

Sample Index: 9

Sample Type: QC

Concentration: 50.0 ng/mL

Calculated Conc: 129. ng/mL

Acq. Date: 11/5/2009

Acq. Time: 4:52:56 PM

Modified: No

Proc. Algorithm: Analyst Classic

Bunching Factor: 1

Noise Threshold: 15.73 cps

Area Threshold: 78.65 cps

Smooths: 1

Sep. Width: 0.20

Sep. Height: 0.01

Exp. Peak Ratio: 5.00

Exp. Adj. Ratio: 4.00

Exp. Val. Ratio: 3.00

RT Window: 30.0 sec

Expected RT: 4.06 min

Use Relative RT: No

Int. Type: Base To Base


Retention Time: 4.00 min

Area: 343599 counts

Height: 77313 cps

Start Time: 3.91 min

End Time: 4.31 min



Sample Name: "L19348-3 Spk C" Sample ID: "Deer #7 3.5 yr male-muscle Spike C, 50 ng/mL, DF=10" File: "110409AD.wiff"

Peak Name: "13CPFOS(M+4)(S)" Mass(es): "503.0/60.0 amu"

Comment: "Annotation: "

Sample Index: 9

Sample Type: QC

Concentration: 1.00 ng/mL

Calculated Conc: N/A

Acq. Date: 11/5/2009

Acq. Time: 4:52:56 PM

Modified: No

Proc. Algorithm: Analyst Classic

Bunching Factor: 1

Noise Threshold: 14.21 cps

Area Threshold: 701.07 cps

Smooths: 1

Sep. Width: 0.20

Sep. Height: 0.01

Exp. Peak Ratio: 5.00

Exp. Adj. Ratio: 4.00

Exp. Val. Ratio: 3.00

RT Window: 30.0 sec

Expected RT: 4.06 min

Use Relative RT: No

Int. Type: Base To Base


Retention Time: 4.00 min

Area: 345013 counts

Height: 122004 cps

Start Time: 3.95 min

End Time: 4.20 min



lyst Version: 1.4.2
rinting Time: 8:44:29 AM
rinting Date: Friday, November 06, 2009

MPI Study: L19346
MPI Set No.: 110409AD

Operator: Mark Neeley
Instrument No.: LC/MS/MS #9

Sample Name: "SS33615" Sample ID: "CCV, 1.0 ng/mL" File: "110409AD.wif"
Peak Name: "PFOS" Mass(es): "499.0/80.0 amu"

Comment: "Annotation"

Sample Index: 10
Sample Type: QC
Concentration: 1.00 ng/mL
Calculated Conc: 1.00 ng/mL 1.00e5
Acq. Date: 11/5/2009 9.50e4
Acq. Time: 5:03:50 PM 9.00e4

Modified: No
Proc. Algorithm: Analyst Classic 8.50e4

Bunching Factor: 1
Noise Threshold: 15.73 cps 8.00e4

Area Threshold: 78.65 cps 7.50e4

a. Smooth: 1 7.50e4

b. Width: 0.20 7.00e4

c. Height: 0.01 6.50e4

d. Peak Ratio: 5.00 6.50e4

Exp. Adj. Ratio: 4.00 6.50e4

Exp. Val. Ratio: 1.00 6.00e4

RT Window: 30.0 sec 5.50e4

Expected RT: 4.06 min 5.50e4

Use Relative RT: No 5.00e4

e. Type: Base To Base 4.50e4

Retention Time: 4.00 min 4.00e4

Area: 365129 counts 4.00e4

Height: 107000 cps 3.50e4

Start Time: 3.00 min 3.00e4

End Time: 4.40 min 3.00e4

2.50e4

2.00e4

1.50e4

1.00e4

5000.00

0.00

4.00
Intensity, cps
2 4 6 8
Time, min

Sample Name: "SS33615" Sample ID: "CCV, 1.0 ng/mL" File: "110409AD.wif"

Peak Name: "13CPFOS(M+X)IST" Mass(es): "503.0/60.0 amu"

Comment: "Annotation"

Sample Index: 10
Sample Type: QC

Concentration: 1.00 ng/mL

Calculated Conc: N/A 1.05e5

Acq. Date: 11/5/2009 1.00e5

Acq. Time: 5:03:50 PM 1.00e5

Modified: No 9.50e4

Proc. Algorithm: Analyst Classic 9.00e4

Bunching Factor: 1 8.50e4

Noise Threshold: 140.21 cps 8.50e4

Area Threshold: 701.07 cps 8.00e4

Num. Peaks: 1 8.00e4

Sep. Width: 0.20 7.50e4

Sep. Height: 0.01 7.00e4

Exp. Peak Ratio: 5.00 7.00e4

Exp. Adj. Ratio: 4.00 6.50e4

Exp. Val. Ratio: 3.00 6.50e4

RT Window: 30.0 sec 6.00e4

Expected RT: 4.06 min 6.00e4

Use Relative RT: No 5.50e4

Int. Type: Base To Base 5.00e4

Retention Time: 4.00 min 4.50e4

Area: 333220 counts 4.00e4

Height: 116407 cps 4.00e4

Start Time: 3.94 min 3.50e4

End Time: 4.20 min 3.50e4

3.00e4

2.50e4

2.00e4

1.50e4

1.00e4

5000.00

0.00

4.00
Intensity, cps
2 4 6 8
Time, min