

U.S. EPA Analytical Methods for the Analysis of Perchlorate in Drinking Water

Dan Hautman

Public Meeting and Webcast

September 20, 2012



Early IC Methods

- **Concern over perchlorate in the early 90's prompted development of early Ion Chromatography (IC) methods**
- **California Dept. of Health Services, June 1997**
 - IC Method with Suppressed Conductivity Detection
 - AS5 analytical column, required eluent modifier (p-cyanophenol)
 - Method was approved for CA State monitoring
 - High Total Dissolved Solids (TDSs) can impact analysis



Early IC Methods

- **New column technologies eliminated need for eluent modifier**
- **Dionex Corporation, 1998, Application Note 121**
 - IC Method with Suppressed Conductivity Detection
 - AS11 analytical column, no eluent modifier needed
 - Adopted as approved alternate for CA State monitoring
 - High TDSs are less of an impact but still can affect analysis




Development of EPA 314.0

- **Further advances in analytical column technology**
- **Interest in establishing a standard EPA test method**
- **EPA Method 314.0, revision 1.0, November 1999**
 - IC Method with suppressed conductivity detection
 - Minimum Reporting Level (MRL) of 4.0 ug/L with detection limit (DL) of 0.5 ug/L
 - High capacity AS16 analytical column, no eluent modifier
 - AS16 further reduces TDS effect
 - Includes additional safeguard to monitor and reduce TDS
 - Sample matrix spike evaluation for confirmation
 - Widely adopted as the standard perchlorate method



Superimposed IC Chromatograms Indicating the Influence of High Concentrations of Common Anions on Low Concentration Measurement of Perchlorate at 4.0 ug/L



Perchlorate Public Meeting and Webcast

5



Additional Method Development Goals

- **Reduce MRL to < 1ug/L through application of:**
 - Sample concentration techniques
 - Microbore analytical columns
 - Advanced detection systems, i.e., mass spectrometry (MS)
- **Further increase tolerance for high ionic strength matrices**
- **Enhance measurement selectivity**

Perchlorate Public Meeting and Webcast

6



Method 314.1, revision 1.0 May 2005

- Lowered the MRL to < 0.2 ug/L (DL 0.03 ug/L) using online sample pre-concentration
- Uses matrix diversion to handle high ionic strength matrices (up to 1000 mg/L TDS)
- Added confirmation analysis using a second analytical column

Perchlorate Public Meeting and Webcast

7



Method 314.2, version 1.0 May 2008

- Lowered the MRL to < 0.1 ug/L (DL 0.02 ug/L) using large volume injection
- Uses 2-D chromatography to handle high ionic strength matrices (up to 1000 mg/L TDS)
- 2-D chromatography eliminates the need for separate confirmation analysis

Perchlorate Public Meeting and Webcast

8



Method 331.0, revision 1.0 January 2005

- **Lowered the MRL to < 0.1 ug/L (DL < 0.01 ug/L)**
- **Applied multiple analytical advancements to an LC analysis**
 - A perchlorate selective LC column (AS-21)
 - MS or MS/MS detection for selectivity and sensitivity
 - A custom labeled internal standard (Cl^{18}O_4)




Method 332.0, revision 1.0 March 2005

- **Lowered the MRL to 0.1 ug/L (DL 0.02 ug/L)**
- **Applied multiple analytical advancements in an IC analysis**
 - Suppressed conductivity IC
 - MS or MS/MS selectivity and sensitivity
 - A custom labeled internal standard (Cl^{18}O_4)



331.0 MS/MS Chromatogram with high TDS



Perchlorate Public Meeting and Webcast

11



Method	Expected Reporting Limit (ug/L)	MDL (ug/L)	Cycle Time Minutes	Acceptable Performance in 1000mg/L TDS	Confirmation	Complexity	Relative Cost
314.0	4.0	0.53	11	No	Matrix Spike Assessment	Moderate	\$45 - \$85
314.1	0.14	0.03	25 (50 with confirmation)	Yes	Second Column Analysis	High	\$60 - \$110
314.2	0.06	0.018	37	Yes	2-D	High	\$60 - \$110
331.0	0.022 (MS/MS)	0.005 (MS/MS)	9	Yes	MS/MS or MS	Moderate	\$80 - \$180
	0.056 (SIM)	0.008 (SIM)					
332.0	0.10 (SIM)	0.02 (SIM)	8	Yes	MS	Moderate	\$80 - \$180

Perchlorate Public Meeting and Webcast

12