From: McCord, James [/O=EXCHANGELABS/OU=EXCHANGE ADMINISTRATIVE GROUP

(FYDIBOHF23SPDLT)/CN=RECIPIENTS/CN=MCCORD, JAMES]

Sent: 5/25/2018 12:13:33 PM

To: Detlef R. U. Knappe [/o=ExchangeLabs/ou=Exchange Administrative Group

(FYDIBOHF23SPDLT)/cn=Recipients/cn=user17c3f77b]; Strynar, Mark [/o=ExchangeLabs/ou=Exchange


Administrative Group (FYDIBOHF23SPDLT)/cn=Recipients/cn=5a9910d5b38e471497bd875fd329a20a-Strynar, Mark]

CC: Hopkins, Zachary [zrhopkin@ncsu.edu]; Nadine Kotlarz [/o=ExchangeLabs/ou=Exchange Administrative Group


(FYDIBOHF23SPDLT)/cn=Recipients/cn=userc79d3fb6]; Chuhui Zhang [czhang24@ncsu.edu]

Subject: RE: Package Delivery


When I was doing structural determination on the Chemours Outfall back in May I was able to resolve the two for the PMPA/PFMOPrA you are discussing. See figure below from the manuscript in prep.



For PEPA/PFMOBA I was only able to detect one chromatographic peak (below)



but the MS/MS contains CF3/CF3O (68.99/84.99) and C2F5/C2F5O (118.99/134.99) (below) which is consistent with a mix of PEPA structure from Chemours and the Sun paper.



From: Detlef Knappe [mailto:knappe@ncsu.edu]

Sent: Friday, May 25, 2018 7:52 AM

To: Strynar, Mark < Strynar. Mark@epa.gov>

Cc: Hopkins, Zachary <zrhopkin@ncsu.edu>; Nadine Kotlarz <nkotlar@ncsu.edu>; Chuhui Zhang <czhang24@ncsu.edu>;

McCord, James <mccord.james@epa.gov>

Subject: Re: Package Delivery

Sure do!

So should the strategy be to sum linear and branched peaks if we get distinguishable signals?

On Fri, May 25, 2018, 7:48 AM Strynar, Mark < Strynar.Mark@epa.gov> wrote:

Strynar, Mark has shared a OneDrive for Business file with you. To view it, click the link below.



PEPA and PMPA structures.pptx

Detlef,

Do you recall when Mei's paper was published we had a discussion about PFMOPrA and PFMOBA being drawn incorrectly for the structures in Figure S1. We did not have standards so we had no way to tell. I believe they are one in the same. PFMOPrA is PMPA and PFMOBA is PEPA. In my 2015 paper I had the PMPA drawn right (Figure S5) but the PEPA wrong (Figure S5).

There MAY still be linear and branched isomers. We can do more QTOF work and run some HILIC methods to see if we get two peaks for these masses. However like PFOS and PFOA have isomers I don't think it is important.

See the attached PPT file

Detlef,

From: zrhopkin@ncsu.edu Sent: Friday, May 25, 2018 6:54 AM To: Detlef R. U. Knappe knappe@ncsu.edu Cc: Nadine Kotlarz nkotlar@ncsu.edu ; Chuhui Zhang czhang24@ncsu.edu ; Strynar, Mark Strynar.Mark@epa.gov ; McCord, James mccord.james@epa.gov > Subject: Re: Package Delivery
If I recall correctly the secondary ion was showing up. This suggested the branched form was present? I would need to compare the MS file parameters for the two in order to see what is different. If they are different then I could test the PFMOBA parameters on the PEPA.
Best,
Zack Hopkins
PhD student
Graduate Research Assistant
Mann Hall 319A Office
Civil, Construction, and Environmental Engineering
North Carolina State University
Raleigh, NC 27695
zrhopkin@ncsu.edu
301-518-7697
On May 24, 2018, at 11:56 PM, Detlef Knappe < knappe@ncsu.edu > wrote: But you were able to see PEPA before, correct? And could distinguish from PFMOBA, I think
On Thu, May 24, 2018 at 11:11 PM, Zachary Hopkins < <u>zrhopkin@ncsu.edu</u> > wrote:

My understanding as 0.1% solution		e compounds PMPA and PEPA. We were provided bot	h of these
I have infused bo detect them.	oth of these on the tr	ple quad. However, once I ran a calibration curve I cou	ıld not
Best,			
Zack			
On Thu. May 24.	2018 at 7:25 PM, De	:lef Knappe < <u>knappe@ncsu.edu</u> > wrote:	
Hi all,	,		
Do we have the	ese compounds?		
The state of the s	TO ANTA MATERIA		
** Water-layer and or adjugate. Yes it was been his some, so done they be did to prove	An analabotea.		

I know we have the linear methoxy versions, but I don't know whether Chemours sent the branched isomers.
Best,
Detlef
On Mon, May 7, 2018 at 5:18 PM, Nadine Kotlarz < nkotlar@ncsu.edu> wrote:
Zack,
The standards we recently received from Chemours were at 1000 ng/uL in water. I prepared 10 ng/uL aqueous stocks for each compound. Then I prepared a 1 ng/uL combined stock with:
1. PFO2HxA
2. NVHOS
3. PFO5DoA
4. PFO4DA
5. PFO3OA
6. PEPA
7. PMPA
8. Nafion byproduct 4
in water.
I already had a 1 ng/uL combined stock of GenX, PFMOAA, Nafion byproduct 1 and Nafion byproduct 2 in methanol. We prepared these standards with solids from Chemours a while ago. We decided we will make our calf serum calibration curves by dosing in the methanol and aqueous mixtures separately.
You're welcome to use any of the dilutions you want. They're stored in the cabinet to the bottom left of Mark's work space in the lab with the triple quad.

	Campus Box 7908 Raleigh, NC 27695				
	Phone: 919-515-7628 Fax: 919-515-7908				
	Best,				
	Zack Hopkins				
	PhD student				
	Graduate Research Assistant				
	Mann Hall 319A Office				
	Civil, Construction, and Environmental Engineering				
	North Carolina State University				
	Raleigh, NC 27695				
	zrhopkin@ncsu.edu				
	301-518-7697				
	Best,				
	Zack Hopkins				
***************************************	PhD student				
	Graduate Research Assistant				

Mann Hall 319A Office

Civil, Construction, and Environmental Engineering

North Carolina State University

Raleigh, NC 27695

zrhopkin@ncsu.edu

301-518-7697