

Message

From: Strynar, Mark [/O=EXCHANGELABS/OU=EXCHANGE ADMINISTRATIVE GROUP (FYDIBOHF23SPDLT)/CN=RECIPIENTS/CN=5A9910D5B38E471497BD875FD329A20A-STRYNAR, MARK]
Sent: 2/13/2019 5:14:49 PM
To: 'Leung, Lam-Wing H' [LAM.H.LEUNG-1@chemours.com]
CC: Swank, Adam [Swank.Adam@epa.gov]; McCord, James [mccord.james@epa.gov]
Subject: RE: E1 GC method

Lam,

Hope all is well with you. I see this email exchange some time back but I don't see the attached GC/MS headspace method. Can you resend it to me?

Thanks,
Mark

From: Strynar, Mark
Sent: Tuesday, November 06, 2018 11:25 AM
To: Leung, Lam-Wing H <LAM.H.LEUNG-1@chemours.com>
Cc: Swank, Adam <Swank.Adam@epa.gov>
Subject: RE: E1 GC method

Thanks a bunch Lam. Are you at SETAC?

Mark

From: Leung, Lam-Wing H [mailto:LAM.H.LEUNG-1@chemours.com]
Sent: Tuesday, November 06, 2018 10:57 AM
To: Strynar, Mark <Strynar.Mark@epa.gov>
Cc: Swank, Adam <Swank.Adam@epa.gov>
Subject: RE: E1 GC method

Hi Mark,

Please attached for the headspace GC/MS method we have been using for E1 analysis. Please note that the method uses SIM for detection (and obviously one can also use the SCAN mode for confirmation purpose). As for the GC column, a shorted column can also be used but if there are other "volatiles" present, it might complicate the separation. Hope this is helpful and please let me know if more info is needed.

Best Regards,
Lam

Lam Leung, Ph.D.
Technical Fellow

lam.h.leung-1@chemours.com

302 353 5218 o
302 985 1655 m

The Chemours Company
Experimental Station 402/5323

200 Powder Mill Road
Wilmington, DE 19803

USA



[LinkedIn](#) | [Twitter](#) | [Chemours.com](#)

From: Strynar, Mark [<mailto:Strynar.Mark@epa.gov>]
Sent: Tuesday, November 06, 2018 10:06 AM
To: Leung, Lam-Wing H <LAM.H.LEUNG-1@chemours.com>
Cc: Swank, Adam <Swank.Adam@epa.gov>
Subject: E1 GC method

Hi Lam,

Hope all is well with you. I know we discussed in the past the disappearance of HFPO-DA in DMSO and the conversion to E1. Some of our colleagues wanted use to measure the evolution of E1 as the HFPO-DA is lost. Do you have a GC method for the measurement of E1 you could share with us.

Mark

Dr. Mark J. Strynar
Physical Scientist
US EPA
National Exposure Research Laboratory
919-541-3706
Strynar.mark@epa.gov

This communication is for use by the intended recipient and contains information that may be privileged, confidential or copyrighted under applicable law. If you are not the intended recipient, you are hereby formally notified that any use, copying or distribution of this e-mail, in whole or in part, is strictly prohibited. Please notify the sender by return e-mail and delete this e-mail from your system. Unless explicitly and conspicuously designated as "E-Contract Intended", this e-mail does not constitute a contract offer, a contract amendment, or an acceptance of a contract offer. This e-mail does not constitute a consent to the use of sender's contact information for direct marketing purposes or for transfers of data to third parties.

Francais Italiano Deutsch Portuges Espanol Japanese Chinese Korean

https://www.chemours.com/Chemours_Home/en_US/email_disclaimer.html