

Message

From: Strynar, Mark [/O=EXCHANGELABS/OU=EXCHANGE ADMINISTRATIVE GROUP (FYDIBOHF23SPDLT)/CN=RECIPIENTS/CN=5A9910D5B38E471497BD875FD329A20A-STRYNAR, MARK]
Sent: 11/16/2020 2:00:33 PM
To: Leung, Lam-Wing H [LAM.H.LEUNG-1@chemours.com]
CC: Risen, Amy J [amy.risen@ncdenr.gov]; Delinsky, Amy [amy.delinsky@ncdenr.gov]
Subject: RE: ?RE: ?isomers of PMPA and PEPA and initial efforts

Hi Lam,

For point 1 I have different retention times for the PMPA and the PFCEA-F so it cannot be spectral artifact as the concentration changes. I can tell them apart. I have run the individual standards and am assessing the content of each analyte in each individual standard today.

For point 2. The well samples were collected in August of 2017. They were during a site visit by the NCDEQ staff. I believe the samples were split between the EPA and either the contract lab for Chemours or Region 4. I am sorry but I don't know those details.

From: Leung, Lam-Wing H <LAM.H.LEUNG-1@chemours.com>
Sent: Thursday, November 12, 2020 11:21 AM
To: Strynar, Mark <Strynar.Mark@epa.gov>
Cc: Risen, Amy J <amy.risen@ncdenr.gov>; Delinsky, Amy <amy.delinsky@ncdenr.gov>
Subject: RE: ?RE: ?isomers of PMPA and PEPA and initial efforts

Hi Mark,

Thanks for the message and the detailed info you provided and I am also delighted to learn that we are in agreement with the "unique" MRMs for both PMPA (229→185) and PEPA (279→235). As for your observation of the 10% "linear PMPA" in the well samples, I have the following comments:

1. Based on our limited experiments (see results box I provided earlier and below for the PMPA only standard), PMPA does show a positive response for the 229→85 MRM and at 1ppb, it represents around 15% so it's probably worthwhile to double check this ration at different concentrations. I also appreciate your observation that the retention time of these peak can be shifting as they elute so early.

MRM transitions	Peak Area				
	PEPA	PMPA	PFCEA A	PFCEA F	Mixed
	1PPB Standard	1PPB Standard	1PPB Standard	1PPB Standard	1PPB Standard
279→234.9	1314.47	No peak	No peak	No peak	1181.01
279→85	No peak	No peak	3381.46	No peak	3495.68
229→184.9	No peak	1494.66	No peak	No peak	1347.27
229→85	No peak	218.6	No peak	3248.8	2987.24

2. Furthermore, can you please provide us with the history of the well samples (specifically when they were collected) such that we can check the exact concentration of PMPA and PEPA concentrations reported for these samples and we can then potentially estimate the "% contribution" of the "linear MRM" based on the PMPA concentration.

Thanks again for the detailed information you provided.

Best Regards,
Lam

Lam Leung, Ph.D.

Technical Fellow
Chemours Discovery Hub N1-106-A
201 Discovery Blvd
Newark, DE 19713

lam.h.leung-1@chemours.com

302 773 6581 o

302 985 1655 m

The Chemours Company



[LinkedIn](#) | [Twitter](#) | [Chemours.com](#)

From: Strynar, Mark <Strynar.Mark@epa.gov>
Sent: Tuesday, November 10, 2020 12:04 PM
To: Leung, Lam-Wing H <LAM.H.LEUNG-1@chemours.com>
Cc: Risen, Amy J <amy.risen@ncdenr.gov>; Delinsky, Amy <amy.delinsky@ncdenr.gov>
Subject: ?RE: ?isomers of PMPA and PEPA and initial efforts

External email. Confirm links and attachments before opening.

Lam,

I have made further progress and can get chromatographic separations here are a few slides I am working on for NCDEQ. (see slide 2 for chromatographic separation). I think we are in agreement based on your text below with some modifications.

I have done work on my QTOF and on my LC-MS/MS system. PFECA-A and PFECA-F DON'T lose a CO₂ as the PMPA and PEPA does easily, so these can be used as unique transitions for PMPA and PEPA.

1. PMPA transitions (after source CO₂ loss) are all unique for 185-119; 185-85; 185-69 (see slide 3)
2. The transition for PFECA-F 229-85 is not unique (see slide 3) but I can get two peaks with the latter being PFECA-F.
3. PEPA transitions (after CO₂ loss) 235-135; 235-119 and 235-69 are all unique (see slide 4).
4. The transition for PFECA-A is unique for 279-85 (see slide 4)

I do see evidence of the liner version of PMPA (PFECA-F) as about 10% of the PMPA peak and a couple of the selected well samples LTW-04 and SMW-06B (see slides 5 and 6) I shown here. All so far show evidence of this peak at about 10% of the PMPA peak. I am working on some integrations to share. The RT of the early eluters are really prone to some movement so the RT do drift about a bit.

For the linear versions of PEPA (PFECA-A) in the real well samples for 279-85 I see two peaks but they are very low responders (no slide shown) an neither co-align with PEPA???

Mark

From: Leung, Lam-Wing H <LAM.H.LEUNG-1@chemours.com>

Sent: Tuesday, November 10, 2020 10:42 AM

To: Strynar, Mark <Strynar.Mark@epa.gov>

Subject: RE: ?isomers of PMPA and PEPA and initial efforts

Hi Mark,

Thanks for your message and I apologize for the late response on this. As I mentioned previously, your "initial finding" is consistent with what we observed. I have listed below results from our experiments using our "Table 3" method for 1ppb standard in water (the mixed standard contains 1ppb of each individual analyte). As you can see, we might not be able to separate them chromatographically, we can successfully separate them spectrally. We are in the process of running the same standards with Method 533 as it uses somewhat different conditions and I'll keep you posted on this. I do like to know if you have made any more progress on this recently and perhaps we can discuss further later this week or sometime next week. We have communicated these limited findings to DEQ and I am hoping that we will be able to "settle" this soon. Thanks again.

MRM transitions	Peak Area				
	PEPA	PMPA	PFECA A	PFECA F	Mixed
	1PPB Standard	1PPB Standard	1PPB Standard	1PPB Standard	1PPB Standard
279→234.9	1314.47	No peak	No peak	No peak	1181.01
279→85	No peak	No peak	3381.46	No peak	3495.68
229→184.9	No peak	1494.66	No peak	No peak	1347.27
229→85	No peak	218.6	No peak	3248.8	2987.24

Best regards,
Lam

Lam Leung, Ph.D.

Technical Fellow

Chemours Discovery Hub N1-106-A

201 Discovery Blvd

Newark, DE 19713

lam.h.leung-1@chemours.com

302 773 6581 o

302 985 1655 m

The Chemours Company



[LinkedIn](#) | [Twitter](#) | [Chemours.com](#)

From: Strynar, Mark <Strynar.Mark@epa.gov>
Sent: Monday, October 26, 2020 10:09 AM
To: Leung, Lam-Wing H <LAM.H.LEUNG-1@chemours.com>
Subject: ?isomers of PMPA and PEPA and initial efforts

External email. Confirm links and attachments before opening.

Hi Lam,

I have done some initial work on my QTOF for looking at PMPA vs PFECA-F (aka PFMOPrA) and for PEPA vs. PFECA-A (aka PFMOBA).

I still need to do work on my QQQ, however on the QTOF here are my initial findings and I wanted to see if they support what you are seeing. I ran each compound in MS mode only then I ran MS/MS at 10, 20 and 40 volts looking at the precursor ion for each M-H, and the CO₂ loss mass. All of these were prepared and run at approximately the same concentrations. The PEPA and PMPA came from the 0.1% in water solutions you gave to me some time back. The PFECA-A and PFECA-F came from neat preparations I made up last week.

1. PMPA (m/z 228.9) readily decarboxylates in the gas phase of the source to give the ion (m/z 184.9). The ratio of the M-CO₂ to the M-H is about 10:1 in the source alone.
2. PEPA (m/z 278.9) likewise readily decarboxylates like PMPA to give the ion (m/z 234.9). The ratio is similar at 10:1 for the M-CO₂ to the M-H is about 10:1 also.
3. The in-source loss of the CO₂ is consistent with what I see with the HFPO-DA (m/z 328.9) and the corresponding fragment (m/z 284.9) so this makes sense I see this for PMPA and PEPA.
4. PFECA-F and PFECA-A really don't seem to lose the CO₂ readily in the source or in MS/MS mode. The M-CO₂ peak for both is really non-existent. When I isolate either and change fragmentor voltage's the main ion is the CF₃O (m/z 84.9) for both. I will have to look for any other fragments more closely.

I will have more later but these are my first impressions.

Mark

Dr. Mark J. Strynar
Physical Scientist
US EPA
National Exposure Research Laboratory
919-541-3706
Strynar.mark@epa.gov

<="" style="">

This communication is for use by the intended recipient and contains information that may be privileged, confidential or copyrighted under applicable law. If you are not the intended recipient, you are hereby formally notified that any use, copying or distribution of this e-mail, in whole or in part, is strictly prohibited. Please notify the sender by return e-mail and delete this e-mail from your system. Unless explicitly and conspicuously designated as "E-Contract Intended", this e-mail does not constitute a contract offer, a contract amendment, or an acceptance of a contract offer. This e-mail does not constitute a consent to the use of sender's contact information for direct marketing purposes or for transfers of data to third parties.

<https://www.chemours.com/en/email-disclaimer>

<="" style="">

This communication is for use by the intended recipient and contains information that may be privileged, confidential or copyrighted under applicable law. If you are not the intended recipient, you are hereby formally notified that any use, copying or distribution of this e-mail, in whole or in part, is strictly prohibited. Please notify the sender by return e-mail and delete this e-mail from your system. Unless explicitly and conspicuously designated as "E-Contract Intended", this e-mail does not constitute a contract offer, a contract amendment, or an acceptance of a contract offer. This e-mail does not constitute a consent to the use of sender's contact information for direct marketing purposes or for transfers of data to third parties.

<https://www.chemours.com/en/email-disclaimer>