



# Per- and Poly-fluorinated compounds Analysis: Case Study Cape Fear Drainage Basin

Mark Strynar USEPA/ORD/NERL/EMMD

Andy Lindstrom, Seth Newton (EPA)  
James McCord, Johnsie Lang (ORISE)




R4 State Labs Presentation  
October 17<sup>th</sup>, 2017



# Presentation Content

- Background
- Timeline
- Methods
- Nakayama et al., 2007 findings
- Strynar et al., 2015 findings
- Sun et al., 2016 findings
- Current efforts and findings

### PFAS Generally Not Removed During Conventional Drinking Water Treatment



Similar for PFOS, PFHxA and PFHxS

Rahman et al., (2014) *Water Research*, 50:318-340


Work approached for detection not quantification

## Timeline: US EPA/NERL Water PFAS Efforts


- 2007 Nakayama et al., PFAS in Cape Fear Watershed
- 2007-2008 Decatur, AL discovery of PFAS in bio-solids, soil, surface water, groundwater
- 2009 (January) Provisional Health Advisory (PHA) for PFOA and PFOS established 400 and 200 ng/L
- 2011-2012 Sampling of Cape Fear River for determination of replacement chemistries
- 2012 (November) Presentation at SETAC Long Beach CA – Identification of novel polyfluorinated compounds in natural waters using accurate mass TOFMS (GenX)
- 2013-2014 Collaborative study with NCSU with sampling of Cape Fear River and three municipalities (Pittsboro, Fayetteville, Wilmington, NC) for legacy and PFECA determination
- 2014 (November) Presentation at SETAC Vancouver, BC - Determination of perfluoroalkyl ether carboxylic acids (PFECA) and sulfonic acids (PFESAs) in North Carolina surface water using high resolution mass spectrometry (GenX and eleven other novel chemicals)
- 2015 (August) Publication of Strynar et al., *"Identification of Novel Perfluoroalkyl Ether Carboxylic Acids (PFECA) and Sulfonic Acids (PFESAs) in Natural Waters Using Accurate Mass Time-of-Flight Mass Spectrometry (TOFMS)"* ES&T
- 2016 (May) OW Health Advisory for PFOS and PFOA set at 70 ng/L
- 2016 (November) Publication of Sun et al., *"Legacy and Emerging Perfluoroalkyl Substances Are Important Drinking Water Contaminants in the Cape Fear River Watershed of North Carolina"* ES&T Letters
- 2017 (June) Reporting of PFECA in CFR and Wilmington drinking water – Star News Online

## NERL/EMMD Water Analysis Method

Surface Water  
Ground Water  
WWTP Effluent  
Drinking Water



## Some Typical PFAS



**Analytes (n=24 native; n=19 isotope labeled)**

Wellington PFAC-24PAR (2 ng/uL) \$250/1.2 mL

	Analyte Name	Acronym	CAS Number
N-1	Perfluorotetradecanoic acid	PFTreA	376-06-7
N-2	Perfluorotridecanoic acid	PFTriA	72629-94-8
N-3	Perfluorododecanoic acid	PFDoA	307-55-1
N-4	Perfluoroundecanoic acid	PFUnA	2058-94-8
N-5	Perfluorodecanoic acid	PFDA	335-76-2
N-6	Perfluorononanoic acid	PFNA	375-95-1
N-7	Perfluoroctanoic acid	PFOA	335-67-1
N-8	Perfluoroheptanoic acid	PFHpA	375-85-9
N-9	Perfluorohexanoic acid	PFHxA	307-24-4
N-10	Perfluoropentanoic acid	PPPeA	2706-90-3
N-11	Perfluorobutanoic acid	PFBA	375-22-4
N-12	Perfluorodecanesulfonate	PFDS	335-77-3
N-13	Perfluorononanesulfonate	PFNS	68259-12-1
N-14	Perfluoroctanesulfonate	PFOS	1763-23-1
N-15	Perfluoroheptanesulfonate	PFHsS	375-92-8
N-16	Perfluorohexanesulfonate	PFHxS	355-46-4
N-17	Perfluoropentansulfonate	PPPeS	2706-91-4
N-18	Perfluorobutanesulfonate	PFBS	375-73-5
N-19	Perfluoroctanesulfonamide	PFOSA	754-91-6
N-20	Fluorotelomer sulfonate 8:2	FtS 8:2	39108-34-4
N-21	Fluorotelomer sulfonate 6:2	FtS 6:2	27619-97-2
N-22	Fluorotelomer sulfonate 4:2	FtS 4:2	NA
N-23	N-ethyl-N-((heptadecafluoroctyl)sulfonyl)glycine	NEtFOSAA	2991-50-6
N-24	N-(Heptadecafluoroctylsulfonyl)-N-methyl/glycine	NMeFOSAA	2355-31-9

**Analytes (n=24 native; n=19 isotope labeled)**

Wellington MPFAC-24ES (1 ng/uL) \$700/1.2 mL

	Analyte Name	Acronym
1-1	PERFLUOROALKYL CARBOXYLATES (PFCAs)	
1-2	Perfluoro-n-[1,2, <sup>13</sup> C <sub>2</sub> ]tetradecanoic acid	M2PFTeDA
1-3	Perfluoro-n-[1,2, <sup>13</sup> C]dodecanoic acid	M2PFDoA
1-4	Perfluoro-n-[1,2,3,4,5,6,7, <sup>13</sup> C]undecanoic acid	M7PFUdA
1-5	Perfluoro-n-[1,2,3,4,5,6, <sup>13</sup> C]decanoic acid	M6PFDA
1-6	Perfluoro-n-[ <sup>13</sup> C]nonanoic acid	M9FNA
1-7	Perfluoro-n-[ <sup>13</sup> C]octanoic acid	M8PFOA
1-8	Perfluoro-n-[1,2,3,4- <sup>13</sup> C]heptanoic acid	M4PFHpa
1-9	Perfluoro-n-[1,2,3,4,6- <sup>13</sup> C]hexanoic acid	M5PFHxA
1-10	Perfluoro-n-[1,2,3,4- <sup>13</sup> C]pentanoic acid	M5PFPeA
1-11	Perfluoro-n-[1,2,3,4- <sup>13</sup> C]butanoic acid	MPFBA
1-12	PERFLUOROALKYLSULFONATES (PFASs)	
1-13	Sodium perfluoro-[ <sup>13</sup> C]octanesulfonate	M8PFOS
1-14	Sodium perfluoro-1-[1,2,3- <sup>13</sup> C]hexanesulfonate	M3PFHxS
1-15	Sodium perfluoro-1-(2,3,4- <sup>13</sup> C)butanesulfonate	M3PFBs
1-16	PERFLUOROOCTANESULFONAMIDES (FOSAs)	
1-17	Perfluoro-1-[ <sup>13</sup> C]octanesulfonamide	M8FOSA-i
1-18	TELOMER SULFONATES	
1-19	Sodium 1H,1H,2H,2H-perfluoro-1-[1,2- <sup>13</sup> C]decane sulfonate (8:2)	M2-8:2FTS
1-20	Sodium 1H,1H,2H,2H-perfluoro-1-[1,2- <sup>13</sup> C]octane sulfonate (6:2)	M2-6:2FTS
1-21	Sodium 1H,1H,2H,2H-perfluoro-1-[1,2- <sup>13</sup> C]hexane sulfonate (4:2)	M2-4:2FTS
1-22	PERFLUOROOCTANESULFONAMIDOACETIC ACIDS	
1-23	N-ethyl-d5-perfluoro-1-octanesulfonamidoacetic acid	d5-N-EtFOSAA
1-24	N-methyl-d3-perfluoro-1-octanesulfonamidoacetic acid	d3-N-MeFOSAA


## Quality Assurance – Quality Control

- Trip spike(s) (ex. 50 ng/L) MOST PEOPLE DO NOT DO THIS
- Trip blank (DI)
- Replication (>10%)
- Occasionally standard addition
- Extracted calibration curves




## Method Summary – Sample Processing


Nakayama et al., 2009; Strynar et al., 2015; Sun et al., 2016; Newton et al., 2017



## Positive Displacement Sample Concentration



CHRATEC Sep-Pak Concentrator  
SPC 10-C (10 ml/min)



Waters Oasis WAX-plus

# Sample Elution



Note: Rack modified to fit 15 mL BD Falcon tubes

Wash:  
4 mL pH 4.0 buffer  
(Acetate/Acetic Acid water)

4 mL methanol

Elute:  
4 mL NH<sub>4</sub>OH/MeOH

Evaporate:  
~15-20 minutes 40°C N<sub>2</sub>  
TurboVap concentrator

~ 0.5 – 1.0 mL

Sample prep for analysis:  
Prepare sample to match  
starting conditions of UPLC  
gradient (75:25 aqu:MeOH)

## Samples That Exceed Calibration Curve

- Estimate concentration based on exceedance analysis
- Dilute sample (analyte and IS presence) with DI lab water (IS added)  
(ex 1:1, 2:1, 5:1 10:1) aiming for midpoint of calibration curve.

Example: Standard curve 10 – 1000 ng/L PFOA  
Sample estimated to be 1600 ng/L


- Dilute sample into new HDPE bottle:  
1 part (250 mL prepared water sample with IS spike and analytes)  
1 part(250 mL) DI water with IS spike at same level as real sample  
Re-extract, elute and analyze.
- Calibrate against existing calibration curve.

# Source of Analytical Standards

Stable Isotope Labeled (<sup>13</sup>C, <sup>18</sup>O, D):

Wellington Labs

Cambridge Isotopes Labs



Native:

Wellington Labs (Mixtures – 1-2 ng/uL; Individual 50 ng/uL)

Sigma-Aldrich

Synquest Labs

Oakwood Products

Manufacturer (industrial mixtures 3M, DuPont,  
Mason Chemicals)



### Method Summary – Water Sample Collection

SS Kemmerer sampler



Lab-made dip sampler



Nalgene 1L HDPE bottles: PFAS free; Prefer wide mouth bottles

## Water Sampling Cape Fear River, NC Nakayama et al., 2007




TABLE 3. Measured Concentrations at the Eleven Sites with the Highest total Concentrations of PFCs in the Cape Fear River Basin\* (See Figure 1 for locations)

site	river	C12 (ng/L)	C11 (ng/L)	C10 (ng/L)	C9 (ng/L)	C8 (ng/L)	C7 (ng/L)	C6 (ng/L)	PFOS (ng/L)	PFHS (ng/L)	PFBs (ng/L)	Total (ng/L)
1	Haw River	9.86	52.1	126	178	287	116	21.7	127	8.43	9.47	842
2	Haw River	3.20	28.7	112	157	200	65.8	12.5	33.4	7.87	2.61	826
3	Haw River	3.29	27.5	103	157	191	69.2	13.7	36.4	9.48	3.04	809
4	Haw River	1.98	20.0	85.2	151	201	66.0	13.2	31.8	7.49	2.88	574
5	Inletury to Cape Fear	2.26	15.0	15.6	15.2	58.6	22.9	23.0	30.5	3.71	ND	531
6	Haw River	1.18	8.87	21.0	72.1	152	65.5	15.5	31.2	7.73	ND	576
7	Cape Fear River	< LOQ	5.34	13.2	34.9	71.3	24.0	6.84	10.7	5.58	ND	227
8	Cape Fear River	1.14	6.38	17.2	35.7	71.8	25.9	6.85	56.4	4.82	ND	228
9	Cape Fear River	1.23	6.75	17.1	33.9	72.7	23.7	7.05	40.7	4.16	ND	211
10	Cape Fear River	< LOQ	7.55	13.3	31.2	46.8	33.9	4.92	56.2	8.84	2.12	199
11	Little River	< LOQ	< LOQ	2.17	2.24	12.6	3.38	3.23	12.2	26.4	3.20	185

\*italicized values show maximal concentrations of each compound.

{Date?Time?}


Framed talk in terms of defining the problem  
Why this is important at NERL

First Sampling 2011-  
2012

Presented at SETAC  
2012

Follow-up sampling  
2013-2013

Presented at SETAC  
2014



## Legacy PFAS found in Cape Fear Water



Analyte	001	002	003	004	005	006	007	008	008	009
C4	23	502	3761	6	4	0	8	7	5	3
C5	441	5607	43590*	17	9	1	32	46	12	9
PFBS	4	5	3	4	5	2	9	5	6	4
C6	17	90	434	18	12	2	27	16	18	14
C7	37	599	3873	14	17	0	11	20	21	9
PFHS	7	12	10	9	7	4	9	10	9	22
C8	32	39	71	33	25	2	38	36	41	18
C9	13	34	127	7	11	1	6	8	11	5
PFOS	19	27	26	17	23	0	0	16	18	14
C10	10	17	12	11	0	3	3	8	10	5

items in red exceed the standard curve high end of 500 ng/l; 10x diluted and re-analyzed; \* still exceed curve and are estimated

## Negative Mass Defect: Per- and Poly- fluorinated Compounds

Formula	Mass	C	H	N	O	P	S	Cl	F
	12.000000								
$\text{C}_8\text{H}_{18}$	114.1408506	8	18						
$\text{C}_8\text{H}_{16}\text{O}_2$	144.1150298	8	16		2				
$\text{C}_8\text{HO}_2\text{F}_{15}$	413.9737023	8	1		2				15
$\text{C}_8\text{HO}_3\text{F}_{15}$	429.9686169	8	1		3				15
$\text{C}_8\text{HO}_3\text{SF}_{17}$	499.937494	8	1		3		1		17

Octane    114.1408506

Octanoic acid    144.1150298

PFOA    413.9737023


Perfluoro-4-ether-octanoic acid    429.9686169

PFOS    499.937494


Mass fragments of per and poly fluorinated compounds conserve negative mass defect

**Identification of Novel Perfluoroalkyl Ether Carboxylic Acids (PFECAs) and Sulfonic Acids (PFESAs) in Natural Waters Using Accurate Mass Time-of-Flight Mass Spectrometry (TOFMS)**


Mark Strynar,<sup>1,2</sup> Sonia Dagnino,<sup>1,3</sup> Rebecca McMahan,<sup>1,4</sup> Shuang Liang,<sup>1,2\*</sup> Andrew Lindstrom,<sup>1</sup> Erik Andersen,<sup>5</sup> Latty McMillan,<sup>6</sup> Michael Thurman,<sup>5</sup> Imma Ferrier,<sup>6</sup> and Carol Ball<sup>3</sup>


**Table 1.** Accurate Mass of Polyfluorinated Compounds and In-Source Artifacts Found in Extracted Water Samples


number	formula	CAS no.	name	[M] <sup>+</sup> <sub>n</sub>	[M + H] <sup>+</sup> <sub>n</sub>	[M + 2H] <sup>2+</sup> <sub>n</sub>	[2M + H] <sup>+</sup> <sub>n</sub>
<b>Monomer PFECAs</b>							
1	C <sub>6</sub> H <sub>12</sub> O <sub>3</sub>			179.9446	178.9773	380.9438	383.9643
2	C <sub>7</sub> H <sub>14</sub> O <sub>3</sub>			219.9813	218.9740	438.9733	458.9553
3	C <sub>8</sub> H <sub>16</sub> O <sub>3</sub>	863690-69-5		279.9762	278.9709	580.9510	558.9491
4	C <sub>9</sub> H <sub>18</sub> O <sub>3</sub>	13232-13-6	unidentified 2-methyl-3-oxahexanoic acid	329.9730	328.9677	680.9247	658.9427
5	C <sub>10</sub> H <sub>20</sub> O <sub>3</sub>			379.9713	378.9645	780.9183	758.9363
6	C <sub>11</sub> H <sub>22</sub> O <sub>3</sub>			438.9885	428.9813	880.9118	858.9293
<b>Polyether PFECAs</b>							
7	C <sub>10</sub> H <sub>18</sub> O <sub>3</sub>	38485-31-8	perfluoro-3,5,7-tri-1-pentenaoctadecanoic acid	445.9513	442.9447	908.8776	886.8957
8	C <sub>11</sub> H <sub>20</sub> O <sub>3</sub>	38487-90-3	perfluoro-3,5,7,9-tetraoctadecanoic acid	477.9588	476.9525	976.8843	754.9113
9	C <sub>12</sub> H <sub>22</sub> O <sub>3</sub>	39492-89-2	perfluoro-3,5,7,9,11-pentaoctadecanoic acid	511.9881	510.9808	1041.9108	822.9289
10	C <sub>13</sub> H <sub>24</sub> O <sub>3</sub>	39493-98-4	perfluoro-3,5,7,9-decahexanoic acid	545.9764	542.9991	1122.9274	990.9435
<b>PFESAs</b>							
11	C <sub>10</sub> H <sub>18</sub> O <sub>5</sub>	66758-36-3 <sup>b</sup>		445.9377	442.9264		
12	C <sub>11</sub> H <sub>20</sub> O <sub>5</sub>			463.9369	462.9326		

**Example Structures**


## m/z Extracted Ion Chromatograph (EIC) Suspected polyfluorinated compound homologous series



## Homologous series NOT offset by CF<sub>2</sub> BUT CF<sub>2</sub>O



**Legacy and Emerging Perfluoroalkyl Substances Are Important  
Drinking Water Contaminants in the Cape Fear River Watershed of  
North Carolina**

Mei Sun,<sup>\*§‡</sup> Elisa Arevalo,<sup>†</sup> Mark Strynar,<sup>§</sup> Andrew Lindstrom,<sup>§</sup> Michael Richardson,<sup>‡</sup> Ben Kearns,<sup>‡</sup> Adam Pickett,<sup>‡</sup> Chris Smith,<sup>#</sup> and Detlef R. U. Knappe<sup>§</sup>

June 7, 2017 story on in Wilmington Star Online News basis


## Toxin taints CFPUA drinking water



### MOST POPULAR

- 1 [Charlotte City teacher, manager of budget - compromised, suspended](#)  
By [Lori L. Johnson](#)
- 2 [North Carolina book, written by local author, goes viral](#)  
By [Troy D. Williams](#)
- 3 [Nags Head couple takes charges](#)  
By [Mike S. Smith](#)
- 4 [Pamlico Seafood, 200000+ likes](#)  
By [Mike S. Smith](#)  
[600,000+ fans](#)

### OUR PICKS



Cape Fear River  
Fayetteville to Wilmington, NC



▲ 1000 CAPTION

A 2000-acre plot of land is up for auction in Fayetteville, NC. The site, home to several ponds, new or older houses, and a school, is about 100 miles upstream from Wilmington. (DID YOU GET ALL THE FACTS? BY LISA WILSON)

By [Lisa Wilson](#) | [lisa.wilson@starnews.com](#)


Published: 10/17/2017 10:17 AM EDT

Updated: 10/17/2017 10:20 AM EDT


Printed: 10/17/2017 10:20 AM EDT

Comments: 0

Share:




Utility can't filter out chemical produced upriver



**Figure 1.** Occurrence of PFASs at drinking water intakes in the CFR watershed. Concentrations represent averages of samples collected between June and December 2013. Individual samples with concentrations below the quantitation limits (QLs) were considered as 0 when calculating averages, and average concentrations below the QLs were not plotted.


Sun et al., 2017

## PFAS removal at Sweeney WTP







Sun et al., 2017

## Other PFECAS Sweeney WTP



■ FFPtOPrA  
FFO2HxA ■ PPMOAA  
FFO3OA ■ PFMOPtA  
FFO4DA ■ PFMoBA




## Current Efforts: EPA/NERL/EMMD


- Completed method validation of PFAS (EPA 24 analyte list) in surface water, waste water and DI water
- Working with NC DEQ and EPA R4 for sample analysis CFR contamination GenX, Nafion by-products and related
- Full method validation GenX and 4 other PFECAs in tap water, surface water, waste water.
- Monitoring well and facility process sample analysis (ongoing)

DUPONT FAYETTEVILLE PLANT  
22828 NC HIGHWAY 87 WEST  
FAYETTEVILLE, NC 28306-7332


[http://iaspub.epa.gov/enviro/tsc/get\\_chem\\_info?v\\_registry\\_id=110000559609](http://iaspub.epa.gov/enviro/tsc/get_chem_info?v_registry_id=110000559609)




### Hexafluoropropylene oxide (HFPO)



- an intermediate used in industrial organofluorine chemistry
- a monomer for fluoropolymers.
- colourless gas is the epoxide of hexafluoropropylene
- produced by DuPont and 3M and as a precursor to the lubricant Krytox and related materials
- generated by oxidation of perfluoropropylene, e.g. with oxygen as well as other oxidants




### PFECAs in 5-15-17 Cape Fear River Sample




ED\_005565\_00008324-00031

5-15-17 EPA Analysis present in CFR



6-23-17 EPA Quality Assurance Spiked Samples



# PFECAs Method Validation

McCord et al., 2017 draft

Figure 1. GenX Method Validation Results. Un-spiked (BL), 50 ng/L (LS) and 200 ng/L (HS) spiked in wastewater, surface water and drinking water.









Figure 2. PFeca A Method Validation Results. Un-spiked (BL), 50 ng/L (LS) and 200 ng/L (HS) spiked in wastewater, surface water and drinking water.




## "Nafion Like" Compounds




Nafion Like Compounds –  
Week#1 Chemours outfall



Nafion Like Compounds  
Week#6 Chemours outfall





**Figure 7. Relative change (compared to highest measured value) in PFAS concentration over weeks 1 – 6 for GenX and NTAs at the Chemours outfall and Sweeney Finished Drinking Water. GenX and NTAs in Panels A,D,E, & F show a consistent decreasing profile. The PFESA Byproduct concentrations are variable and do not show a clear trend.**

Questions?

strynar.mark@epa.gov

Deep River, Moncure, NC