

## Data Verification Module (DVM)

The DVM is an internal review process used by the ADQM group to assist with the determination of data usability. The electronic data deliverables received from the laboratory are loaded into the Locus EIM™ database and processed through a series of data quality checks, which are a combination of software (Locus EIM™ database Data Verification Module (DVM)) and manual reviewer evaluations. The data is evaluated against the following data usability checks:

- Field and laboratory blank contamination
- US EPA hold time criteria
- Missing Quality Control (QC) samples
- Matrix spike(MS)/matrix spike duplicate (MSD) recoveries and the relative percent differences (RPDs) between these spikes
- Laboratory control sample(LCS)/control sample duplicate (LCSD) recoveries and the RPD between these spikes
- Surrogate spike recoveries for organic analyses
- RPD between field duplicate sample pairs
- RPD between laboratory replicates for inorganic analyses
- Difference / percent difference between total and dissolved sample pairs.

There are two qualifier fields in EIM:

**Lab Qualifier** is the qualifier assigned by the lab and may not reflect the usability of the data. This qualifier may have many different meanings and can vary between labs and over time within the same lab. Please refer to the laboratory report for a description of the lab qualifiers. As they are lab descriptors they are not to be used when evaluating the data.

**Validation Qualifier** is the 3rd party formal validation qualifier if this was performed. Otherwise this field contains the qualifier resulting from the ADQM DVM review process. This qualifier assesses the usability of the data and may not equal the lab qualifier. The DVM applies the following data evaluation qualifiers to analysis results, as warranted:

Qualifier	Definition
B	Not detected substantially above the level reported in the laboratory or field blanks.
R	Unusable result. Analyte may or may not be present in the sample.
J	Analyte present. Reported value may not be accurate or precise.
UJ	Not detected. Reporting limit may not be accurate or precise.

The **Validation Status Code** field is set to “DVM” if the ADQM DVM process has been performed. If the DVM has not been run, the field will be blank.

If the DVM has been run (**Validation Status Code** equals “DVM”), use the **Validation Qualifier**.

## DVM Narrative Report

The DVM process was executed, but no data qualifiers were applied for this project.

## ANALYTICAL REPORT


Job Number: 280-106447-1

Job Description: FAY-OLD OUTFALL SAMP 2/18

For:

Chemours Company FC, LLC The  
c/o AECOM  
Sabre Building, Suite 300  
4051 Ogletown Road  
Newark, DE 19713

Attention: Michael Aucoin



Approved for release.  
Michelle A. Johnston  
Project Manager II  
2/22/2018 8:03 AM

---

Michelle A Johnston, Project Manager II  
4955 Yarrow Street, Arvada, CO, 80002  
(303)736-0110  
michelle.johnston@testamericainc.com  
02/22/2018

cc: Barbara McGraw  
Kelly Rinehimer

The test results in this report relate only to the samples in this report and meet all requirements of NELAC, with any exceptions noted. Pursuant to NELAP, this report shall not be reproduced except in full, without the written approval of the laboratory. All questions regarding this report should be directed to the TestAmerica Denver Project Manager.

The Lab Certification ID# is 4025.

Reporting limits are adjusted for sample size used, dilutions and moisture content if applicable.

**TestAmerica Laboratories, Inc.**

TestAmerica Denver 4955 Yarrow Street, Arvada, CO 80002  
Tel (303) 736-0100 Fax (303) 431-7171 [www.testamericainc.com](http://www.testamericainc.com)

# Table of Contents

Cover Title Page . . . . .	1
Data Summaries . . . . .	4
Definitions . . . . .	4
Case Narrative . . . . .	5
Detection Summary . . . . .	8
Client Sample Results . . . . .	13
Default Detection Limits . . . . .	25
Surrogate Summary . . . . .	26
QC Sample Results . . . . .	27
QC Association . . . . .	29
Chronicle . . . . .	30
Certification Summary . . . . .	34
Method Summary . . . . .	35
Sample Summary . . . . .	36
Manual Integration Summary . . . . .	37
Reagent Traceability . . . . .	38
COAs . . . . .	42
Organic Sample Data . . . . .	57
LCMS . . . . .	57
8321A_HFPO_Du . . . . .	57
8321A_HFPO_Du QC Summary . . . . .	58
8321A_HFPO_Du Sample Data . . . . .	65
Standards Data . . . . .	113
8321A_HFPO_Du ICAL Data . . . . .	113
8321A_HFPO_Du CCAL Data . . . . .	138
Raw QC Data . . . . .	156

# Table of Contents

8321A_HFPO_Du Tune Data .....	156
8321A_HFPO_Du Blank Data .....	161
8321A_HFPO_Du LCS/LCSD Data .....	169
8321A_HFPO_Du MS/MSD Data .....	186
8321A_HFPO_Du Duplicate/Triplicate Data .....	190
8321A_HFPO_Du Run Logs .....	194
8321A_HFPO_Du Prep Data .....	196
Shipping and Receiving Documents .....	201
Client Chain of Custody .....	202
Sample Receipt Checklist .....	204

## Definitions/Glossary

Client: Chemours Company FC, LLC The  
Project/Site: FAY-OLD OUTFALL SAMP 2/18

TestAmerica Job ID: 280-106447-1

### Qualifiers

#### LCMS

Qualifier	Qualifier Description
D	Sample results are obtained from a dilution; the surrogate or matrix spike recoveries reported are calculated from diluted samples.
4	MS, MSD: The analyte present in the original sample is greater than 4 times the matrix spike concentration; therefore, control limits are not applicable.

### Glossary

Abbreviation	These commonly used abbreviations may or may not be present in this report.
▫	Listed under the "D" column to designate that the result is reported on a dry weight basis
%R	Percent Recovery
CFL	Contains Free Liquid
CNF	Contains No Free Liquid
DER	Duplicate Error Ratio (normalized absolute difference)
Dil Fac	Dilution Factor
DL	Detection Limit (DoD/DOE)
DL, RA, RE, IN	Indicates a Dilution, Re-analysis, Re-extraction, or additional Initial metals/anion analysis of the sample
DLC	Decision Level Concentration (Radiochemistry)
EDL	Estimated Detection Limit (Dioxin)
LOD	Limit of Detection (DoD/DOE)
LOQ	Limit of Quantitation (DoD/DOE)
MDA	Minimum Detectable Activity (Radiochemistry)
MDC	Minimum Detectable Concentration (Radiochemistry)
MDL	Method Detection Limit
ML	Minimum Level (Dioxin)
NC	Not Calculated
ND	Not Detected at the reporting limit (or MDL or EDL if shown)
PQL	Practical Quantitation Limit
QC	Quality Control
RER	Relative Error Ratio (Radiochemistry)
RL	Reporting Limit or Requested Limit (Radiochemistry)
RPD	Relative Percent Difference, a measure of the relative difference between two points
TEF	Toxicity Equivalent Factor (Dioxin)
TEQ	Toxicity Equivalent Quotient (Dioxin)

**CASE NARRATIVE**  
**Client: The Chemours Company FC, LLC**  
**Project: FAY-OLD OUTFALL SAMP 2/18**  
**Report Number: 280-106447-1**

With the exceptions noted as flags or footnotes, standard analytical protocols were followed in the analysis of the samples and no problems were encountered or anomalies observed. In addition all laboratory quality control samples were within established control limits, with any exceptions noted below. Each sample was analyzed to achieve the lowest possible reporting limit within the constraints of the method. In some cases, due to interference or analytes present at high concentrations, samples were diluted. For diluted samples, the reporting limits are adjusted relative to the dilution required.

Calculations are performed before rounding to avoid round-off errors in calculated results.

All holding times were met and proper preservation noted for the methods performed on these samples, unless otherwise detailed in the individual sections below.

For samples requiring analysis at a dilution, the dilution factor has been multiplied by the Method Detection Limit (MDL) for each analyte and evaluated versus the project-specific reporting limit (PSRL). If the obtained value is below the PSRL, then the PSRL is preserved as the reporting limit for the diluted result, otherwise, the obtained value becomes the reporting limit. This is done in order to maintain the PSRL to meet project requirements at the request of the client and to report the lowest possible RL for each analyte.

#### **Receipt**

The samples were received on 2/15/2018 9:30 AM; the samples arrived in good condition, properly preserved and, where required, on ice. The temperature of the cooler at receipt was 2.7° C.

#### **Receipt Exceptions**

Per client instructions received on 2/15/2018 the job description was logged as FAY-OLD OUTFALL SAMP 2/18.

Per client request for expedited turn around time and laboratory operations approval the requested analyses were logged on a 5 business day turn around time.

No other anomalies were observed during sample receipt.

#### **Standards**

Analytical standards were prepared using the acid form of the compound Perfluoro(2-propoxypropanoic) acid (HFPO-DA).

The surrogate compound, <sup>13</sup>C<sub>3</sub> HFPO-DA was introduced at the extraction step and was used as an internal standard for quantitation of HFPO-DA. The concentration of the surrogate spike is 0.2ug/L in water samples or 50ug/kg in soil samples.

#### **Sample Extraction and Analysis**

The samples presented in this report were extracted for the target analyte by TestAmerica Denver's SOP DV-OP-0019, Rev. 8 and analyzed for the target analyte by TestAmerica Denver's SOP DV-LC-0012, Rev. 14, with the exceptions of the items indicated in the DuPont QAS. Sample FAY-VES-OLDOUTFALL-A (280-106447-2) was chosen to be analyzed as a duplicate and also to be spiked with the target analyte.

For water samples a 250mL aliquot of each sample is extracted using solid phase extraction technique with methanol conditioned Weak Anion Exchange cartridges. Each sample is spiked with the internal standard/surrogate, prior to extraction. After the sample is passed through the cartridge, the analytes are eluted with 2%Formic Acid, 6mLs of HPLC grade MeOH and then with 4mL of 10% ammonium hydroxide in methanol. The final volume is brought to 5mL using reagent water and the extract is analyzed by LC/MS/MS.

The target analyte is separated from other components on a high-performance liquid chromatography (HPLC) C18 column with a mobile phase mixture of water containing 0.1% ammonium acetate and methanol. The mass spectrometer detector is operated in the electrospray (ESI) negative ion mode. The instrument is calibrated at 7 concentration levels (0.2, 0.5, 1.0, 2.0, 5.0, 10 and 20ug/L). The target analyte is detected as the perfluoro(2-propoxypropanoic) acid with the parent ion of 328.8 amu. The daughter ions used for analysis by LC/MS/MS are at 284.8 amu. The ratio of the peak areas to the two ions must be ±20% of the ion ratios in the mid-point ICAL for qualitative identification. Sample results are quantitated using the internal standard dilution.

#### **Tuning and Calibration**

The instrument is tuned with a solution of the target analyte such that mass assignments are within ±0.5 amu of the daughter ions. The instrument is calibrated with seven concentration levels from 0.2ug/L to 20ug/L. Linear regression ( $y=ax+b$ ) or quadratic functions ( $y=ax+cx^2+b$ ) are used with a correlation coefficient or coefficient of determination ≥0.990.

Following initial calibration (ICAL), an initial calibration blank (ICB) is tested, which consists of methanol spiked with the surrogate. The result for the target analyte must be less than one half the reporting limit (RL) to proceed.

Next an initial calibration verification (ICV) standard is tested. This is a mid-level concentration standard from a different vendor from the ICAL standard. If a different vendor is not available then, a different lot number from the same vendor is used. The ICV must be within 80-120% of the true value.

The quantitation limit verification standard is a standard from the same source as the ICAL tested run at the RL level to determine accuracy near the detection limit. This recovery must be within 70-130%.

Continuing calibration verification (CCV) standards are tested every 10 injections and are from the same source as the ICAL and are at mid-level concentration. The recovery of the CCVs must be 70-130% or recalibration is necessary.

#### **Method QC Samples**

The Method Blank is processed reagent water spiked with internal standard and prepared with each batch of 20 samples of the same matrix. All samples in the batch are processed at the same time and with the same reagents. The method blank must be less than the LOD or associated batch samples must be re-extracted and reanalyzed.

Each batch is prepared with a low- and a mid-level concentration spike Laboratory Control Samples (LCS). The recoveries of these samples must be within 70-130% or associated batch samples must be re-extracted and reanalyzed. If the recovery is biased high and samples are non-detect, results can be reported without re-extraction.

#### **Calculations**

##### **Sample Result Calculation**

For internal standard quantitation,

HFPO-DA Response = Area of HFPO-DA \* 13C3 HFPO-DA concentration / area of 13C3 HFPO-DA

Concentration in waters, ug/L = (Cex Vt)/(Vo)

Where:

Cex = Concentration measured in sample extract from the target analyte response (ng/mL)

Vt = Volume of total extract (mL)

Vo = Volume of water extracted (mL)

##### **2. Percent Recovery Calculation**

Spike Recovery = (SSR-SR)/(SA)x100%

Where:

SSR = Spike sample result

SR = Sample result

SA = Spike added

##### **3. Relative Percent Difference Calculation**

RPD = (SR - DR)/(1/2(SR+DR))x100

Where:

SR = Sample result

DR = Duplicate result

#### **HFPO-DA Analysis Anomalies**

Samples FAY-VES-OLDOUTFALL-A-D (280-106447-1), FAY-VES-OLDOUTFALL-A (280-106447-2), FAY-VES-OLDOUTFALL-B (280-106447-3), FAY-VES-OLDOUTFALL-C (280-106447-4), FAY-VES-OLDOUTFALL-D (280-106447-5), FAY-VES-OLDOUTFALL-E (280-106447-6), FAY-VES-OLDOUTFALLSEEP-A (280-106447-7), FAY-VES-OLDOUTFALLCREEK-A (280-106447-8), FAY-VES-OLDOUTFALLCREEK-A2 (280-106447-9), FAY-VES-OLDOUTFALLCREEK-A3 (280-106447-10), FAY-VES-OLDOUTFALLCREEKWATERBO (280-106447-11) and FAY-VES-FB-021418 (280-106447-12) were analyzed for Perfluorinated Hydrocarbons in accordance with DV-LC-0012. The samples were prepared on 02/20/2018 and analyzed on 02/21/2018.

Calibration 9 (STD125) has been included in the raw data, but was not used in the Initial Calibration (ICAL).

Reporting limits have been adjusted accordingly for the initial volumes extracted.

The following sample was decanted prior to preparation: FAY-VES-OLDOUTFALLSEEP-A (280-106447-7).

Each sample is analyzed to achieve the lowest possible reporting limits within the constraints of the method. Due to analytes present above the calibration range, samples FAY-VES-OLDOUTFALL-A-D (280-106447-1), FAY-VES-OLDOUTFALL-A (280-106447-2), FAY-VES-OLDOUTFALL-A (280-106447-2[DU]), FAY-VES-OLDOUTFALL-A (280-106447-2[MS]), FAY-VES-OLDOUTFALL-B (280-106447-3), FAY-VES-OLDOUTFALL-C (280-106447-4), FAY-VES-OLDOUTFALL-D (280-106447-5), FAY-VES-OLDOUTFALL-E (280-106447-6), FAY-VES-OLDOUTFALLCREEK-A (280-106447-8), FAY-VES-OLDOUTFALLCREEK-A2 (280-106447-9), FAY-VES-OLDOUTFALLCREEK-A3 (280-106447-10) and FAY-VES-OLDOUTFALLCREEKWATERBO (280-106447-11) had to be analyzed at dilutions. The reporting limits have been adjusted relative to the dilutions required.

The MS and Sample Duplicate associated with prep batch 280-405473 were performed on sample FAY-VES-OLDOUTFALL-A (280-106447-2). The MS spike compound recovery and surrogate recovery could not be reliably calculated for HFPO-DA and 13C3 HFPO-DA because the sample concentrations were greater than four times the spike amounts. The acceptable low-level LCS and mid-level LCS/LCSD analysis data indicated that the analytical system was operating within control; therefore, corrective action was



deemed unnecessary.

No other analytical or quality issues were noted, other than those described in the Definitions/Glossary page.

## Fluorochemical Characterization of Water Samples Analytical Results

Chemours Sample Identification	TestAmerica Sample Identification	Collection Date/Time	Date Sample Received by TestAmerica	Analysis Date	HFPO-DA# (ug/L**)
FAY-VES-OLDOUTFALL-A-D	280-106447-1	2/14/2018 12:00	2/15/2018	2/21/2018	7.1
FAY-VES-OLDOUTFALL-A	280-106447-2	2/14/2018 12:00	2/15/2018	2/21/2018	7.6
FAY-VES-OLDOUTFALL-B	280-106447-3	2/14/2018 12:45	2/15/2018	2/21/2018	7.3
FAY-VES-OLDOUTFALL-C	280-106447-4	2/14/2018 13:38	2/15/2018	2/21/2018	9.7
FAY-VES-OLDOUTFALL-D	280-106447-5	2/14/2018 14:05	2/15/2018	2/21/2018	13
FAY-VES-OLDOUTFALL-E	280-106447-6	2/14/2018 14:15	2/15/2018	2/21/2018	11
FAY-VES-OLDOUTFALLSEE P-A	280-106447-7	2/14/2018 13:10	2/15/2018	2/21/2018	1.3
FAY-VES-OLDOUTFALLCRE EK-A	280-106447-8	2/14/2018 13:22	2/15/2018	2/21/2018	2.3
FAY-VES-OLDOUTFALLCRE EK-A2	280-106447-9	2/14/2018 16:30	2/15/2018	2/21/2018	2.5
FAY-VES-OLDOUTFALLCRE EK-A3	280-106447-10	2/14/2018 16:45	2/15/2018	2/21/2018	2.2
FAY-VES-OLDOUTFALLCRE EKWATERBO	280-106447-11	2/14/2018 17:00	2/15/2018	2/21/2018	1.8
FAY-VES-FB-021418	280-106447-12	2/14/2018 7:00	2/15/2018	2/21/2018	<0.010

# HFPO-DA – hexafluoropropylene oxide dimer acid, analyzed by method DV-LC-0012, Revision 14.

< = less than the stated value

\*\* ug/L – micrograms/liter (parts per billion)

\*\*\* The MS spike compound recoveries and surrogate recoveries could not be reliably calculated because the sample concentration was greater than four times the spike amounts.

### DEFINITIONS:

Reporting Limit (RL) for the procedure is approximately 0.010ug/L.

### RESULTS ARE CALCULATED ACCORDING TO THE FOLLOWING CRITERIA:

For samples analyzed in duplicate:

If the sample and laboratory duplicate are greater than 5X RL, the relative percent difference (RPD) is less than 20, the average value is reported. If the RPD is greater than 20, the higher value is reported.


If the sample or laboratory duplicate are less than 5X RL, and the absolute difference is less than RL, the average value is reported. If the absolute difference is greater than the RL, the higher value is reported.

### Matrix Spike Recoveries:

Acceptable Range: 70%-130%

TestAmerica Sample ID	Matrix Spike Recoveries
280-106447-2	163%***

**SUBMITTED BY:**



2/22/2018

Michelle A. Johnston, Project Manager

Date

## Executive Summary

Client: Chemours Company FC, LLC The

Job Number: 280-106447-1

### 8321A : HFPO-DA

Lab Sample ID	Client Sample ID	Analyte	Individual Result (ug/L)	Final Result (ug/L)	RL
280-106447-1	FAY-VES-OLDOUTFALL-A-D	HFPO-DA	7.1	7.1	0.044
280-106447-2	FAY-VES-OLDOUTFALL-A	HFPO-DA	7.4	7.6	0.043
280-106447-2 DU	FAY-VES-OLDOUTFALL-A	HFPO-DA	7.8		0.045
280-106447-3	FAY-VES-OLDOUTFALL-B	HFPO-DA	7.3	7.3	0.044
280-106447-4	FAY-VES-OLDOUTFALL-C	HFPO-DA	9.7	9.7	0.047
280-106447-5	FAY-VES-OLDOUTFALL-D	HFPO-DA	13	13	0.045
280-106447-6	FAY-VES-OLDOUTFALL-E	HFPO-DA	11	11	0.045
280-106447-7	FAY-VES-OLDOUTFALLSEEP-A	HFPO-DA	1.3	1.3	0.010
280-106447-8	FAY-VES-OLDOUTFALLCREEK-A	HFPO-DA	2.3	2.3	0.010
280-106447-9	FAY-VES-OLDOUTFALLCREEK-A2	HFPO-DA	2.5	2.5	0.010
280-106447-10	FAY-VES-OLDOUTFALLCREEK-A3	HFPO-DA	2.2	2.2	0.010
280-106447-11	FAY-VES-OLDOUTFALLCREEKWAT ERBO	HFPO-DA	1.8	1.8	0.010
280-106447-12	FAY-VES-FB-021418	HFPO-DA	<0.010	<0.010	0.010

(a) Method 8321A

(b) DUP or REP indicates a laboratory duplicate.

(c) If the sample and laboratory duplicate are both greater than 5X the RL and the relative percent difference (RPD) is less than 20, the average value is reported. If the RPD is greater than 20, the higher of the sample and laboratory duplicate value is reported. If the sample and/or laboratory duplicate are less than 5X the RL, and the absolute difference between the sample and laboratory duplicate is less than the RL, the average value is reported. If the absolute difference is greater than the RL, the higher of the sample and laboratory duplicate value is reported. If either the sample or the duplicate result is greater than or equal to the RL and the other is less than the RL, then the higher of the two is reported.

(d) Moisture Determined by ASTM D2216.

(e) Reporting Limit (RL) = The concentration equivalent to the low calibration standard.

## Detection Summary

Client: Chemours Company FC, LLC The  
Project/Site: FAY-OLD OUTFALL SAMP 2/18

TestAmerica Job ID: 280-106447-1

Client Sample ID: FAY-VES-OLDOUTFALL-A-D

Lab Sample ID: 280-106447-1

Analyte	Result	Qualifier	RL	MDL	Unit	Dil Fac	D	Method	Prep Type
HFPO-DA	7.1		0.044		ug/L	10		8321A	Total/NA

Client Sample ID: FAY-VES-OLDOUTFALL-A

Lab Sample ID: 280-106447-2

Analyte	Result	Qualifier	RL	MDL	Unit	Dil Fac	D	Method	Prep Type
HFPO-DA	7.4		0.043		ug/L	10		8321A	Total/NA

Client Sample ID: FAY-VES-OLDOUTFALL-B

Lab Sample ID: 280-106447-3

Analyte	Result	Qualifier	RL	MDL	Unit	Dil Fac	D	Method	Prep Type
HFPO-DA	7.3		0.044		ug/L	10		8321A	Total/NA

Client Sample ID: FAY-VES-OLDOUTFALL-C

Lab Sample ID: 280-106447-4

Analyte	Result	Qualifier	RL	MDL	Unit	Dil Fac	D	Method	Prep Type
HFPO-DA	9.7		0.047		ug/L	10		8321A	Total/NA

Client Sample ID: FAY-VES-OLDOUTFALL-D

Lab Sample ID: 280-106447-5

Analyte	Result	Qualifier	RL	MDL	Unit	Dil Fac	D	Method	Prep Type
HFPO-DA	13		0.045		ug/L	10		8321A	Total/NA

Client Sample ID: FAY-VES-OLDOUTFALL-E

Lab Sample ID: 280-106447-6

Analyte	Result	Qualifier	RL	MDL	Unit	Dil Fac	D	Method	Prep Type
HFPO-DA	11		0.045		ug/L	10		8321A	Total/NA

Client Sample ID: FAY-VES-OLDOUTFALLSEEP-A

Lab Sample ID: 280-106447-7

Analyte	Result	Qualifier	RL	MDL	Unit	Dil Fac	D	Method	Prep Type
HFPO-DA	1.3		0.010		ug/L	1		8321A	Total/NA

Client Sample ID: FAY-VES-OLDOUTFALLCREEK-A

Lab Sample ID: 280-106447-8

Analyte	Result	Qualifier	RL	MDL	Unit	Dil Fac	D	Method	Prep Type
HFPO-DA	2.3		0.010		ug/L	2		8321A	Total/NA

Client Sample ID: FAY-VES-OLDOUTFALLCREEK-A2

Lab Sample ID: 280-106447-9

Analyte	Result	Qualifier	RL	MDL	Unit	Dil Fac	D	Method	Prep Type
HFPO-DA	2.5		0.010		ug/L	2		8321A	Total/NA

Client Sample ID: FAY-VES-OLDOUTFALLCREEK-A3

Lab Sample ID: 280-106447-10

Analyte	Result	Qualifier	RL	MDL	Unit	Dil Fac	D	Method	Prep Type
HFPO-DA	2.2		0.010		ug/L	2		8321A	Total/NA

Client Sample ID: FAY-VES-OLDOUTFALLCREEKWATERBO

Lab Sample ID: 280-106447-11

—

This Detection Summary does not include radiochemical test results.

TestAmerica Denver

Detection Summary

Client: Chemours Company FC, LLC The  
Project/Site: FAY-OLD OUTFALL SAMP 2/18

TestAmerica Job ID: 280-106447-1

Client Sample ID: FAY-VES-OLDOUTFALLCREEKWATERBO  
(Continued) Lab Sample ID: 280-106447-11

Analyte	Result	Qualifier	RL	MDL	Unit	Dil Fac	D	Method	Prep Type
HFPO-DA	1.8		0.010		ug/L	2		8321A	Total/NA

Client Sample ID: FAY-VES-FB-021418 Lab Sample ID: 280-106447-12

No Detections.

This Detection Summary does not include radiochemical test results.

# Client Sample Results

Client: Chemours Company FC, LLC The  
Project/Site: FAY-OLD OUTFALL SAMP 2/18

TestAmerica Job ID: 280-106447-1

Client Sample ID: FAY-VES-OLDOUTFALL-A-D

Lab Sample ID: 280-106447-1

Date Collected: 02/14/18 12:00

Matrix: Water

Date Received: 02/15/18 09:30

## Method: 8321A - HFPO-DA

Analyte	Result	Qualifier	RL	MDL	Unit	D	Prepared	Analyzed	Dil Fac
HFPO-DA	7.1		0.044		ug/L		02/20/18 10:22	02/21/18 08:50	10
Surrogate	%Recovery	Qualifier	Limits				Prepared	Analyzed	Dil Fac
13C3 HFPO-DA	92	D	50 - 200				02/20/18 10:22	02/21/18 08:50	10

TestAmerica Denver

# Client Sample Results

Client: Chemours Company FC, LLC The  
Project/Site: FAY-OLD OUTFALL SAMP 2/18

TestAmerica Job ID: 280-106447-1

Client Sample ID: FAY-VES-OLDOUTFALL-A

Lab Sample ID: 280-106447-2

Date Collected: 02/14/18 12:00

Matrix: Water

Date Received: 02/15/18 09:30

## Method: 8321A - HFPO-DA

Analyte	Result	Qualifier	RL	MDL	Unit	D	Prepared	Analyzed	Dil Fac
HFPO-DA	7.4		0.043		ug/L		02/20/18 10:22	02/21/18 08:54	10
Surrogate	%Recovery	Qualifier	Limits				Prepared	Analyzed	Dil Fac
13C3 HFPO-DA	89	D	50 - 200				02/20/18 10:22	02/21/18 08:54	10

TestAmerica Denver



# Client Sample Results

Client: Chemours Company FC, LLC The  
Project/Site: FAY-OLD OUTFALL SAMP 2/18

TestAmerica Job ID: 280-106447-1

Client Sample ID: FAY-VES-OLDOUTFALL-B

Lab Sample ID: 280-106447-3

Date Collected: 02/14/18 12:45

Matrix: Water

Date Received: 02/15/18 09:30

## Method: 8321A - HFPO-DA

Analyte	Result	Qualifier	RL	MDL	Unit	D	Prepared	Analyzed	Dil Fac
HFPO-DA	7.3		0.044		ug/L		02/20/18 10:22	02/21/18 09:04	10
Surrogate	%Recovery	Qualifier	Limits				Prepared	Analyzed	Dil Fac
13C3 HFPO-DA	90	D	50 - 200				02/20/18 10:22	02/21/18 09:04	10

TestAmerica Denver

# Client Sample Results

Client: Chemours Company FC, LLC The  
Project/Site: FAY-OLD OUTFALL SAMP 2/18

TestAmerica Job ID: 280-106447-1

Client Sample ID: FAY-VES-OLDOUTFALL-C

Lab Sample ID: 280-106447-4

Date Collected: 02/14/18 13:38

Matrix: Water

Date Received: 02/15/18 09:30

## Method: 8321A - HFPO-DA

Analyte	Result	Qualifier	RL	MDL	Unit	D	Prepared	Analyzed	Dil Fac
HFPO-DA	9.7		0.047		ug/L		02/20/18 10:22	02/21/18 09:07	10
Surrogate	%Recovery	Qualifier	Limits				Prepared	Analyzed	Dil Fac
13C3 HFPO-DA	86	D	50 - 200				02/20/18 10:22	02/21/18 09:07	10

TestAmerica Denver

# Client Sample Results

Client: Chemours Company FC, LLC The  
Project/Site: FAY-OLD OUTFALL SAMP 2/18

TestAmerica Job ID: 280-106447-1

Client Sample ID: FAY-VES-OLDOUTFALL-D

Lab Sample ID: 280-106447-5

Date Collected: 02/14/18 14:05

Matrix: Water

Date Received: 02/15/18 09:30

## Method: 8321A - HFPO-DA

Analyte	Result	Qualifier	RL	MDL	Unit	D	Prepared	Analyzed	Dil Fac
HFPO-DA	13		0.045		ug/L		02/20/18 10:22	02/21/18 09:10	10
Surrogate	%Recovery	Qualifier	Limits				Prepared	Analyzed	Dil Fac
13C3 HFPO-DA	86	D	50 - 200				02/20/18 10:22	02/21/18 09:10	10

TestAmerica Denver

# Client Sample Results

Client: Chemours Company FC, LLC The  
Project/Site: FAY-OLD OUTFALL SAMP 2/18

TestAmerica Job ID: 280-106447-1

Client Sample ID: FAY-VES-OLDOUTFALL-E

Lab Sample ID: 280-106447-6

Date Collected: 02/14/18 14:15

Matrix: Water

Date Received: 02/15/18 09:30

## Method: 8321A - HFPO-DA

Analyte	Result	Qualifier	RL	MDL	Unit	D	Prepared	Analyzed	Dil Fac
HFPO-DA	11		0.045		ug/L		02/20/18 10:22	02/21/18 09:13	10
Surrogate	%Recovery	Qualifier	Limits				Prepared	Analyzed	Dil Fac
13C3 HFPO-DA	91	D	50 - 200				02/20/18 10:22	02/21/18 09:13	10

# Client Sample Results

Client: Chemours Company FC, LLC The  
Project/Site: FAY-OLD OUTFALL SAMP 2/18

TestAmerica Job ID: 280-106447-1

Client Sample ID: FAY-VES-OLDOUTFALLSEEP-A

Lab Sample ID: 280-106447-7

Date Collected: 02/14/18 13:10

Matrix: Water

Date Received: 02/15/18 09:30

## Method: 8321A - HFPO-DA

Analyte	Result	Qualifier	RL	MDL	Unit	D	Prepared	Analyzed	Dil Fac
HFPO-DA	1.3		0.010		ug/L		02/20/18 10:22	02/21/18 08:24	1

Surrogate	%Recovery	Qualifier	Limits	Prepared	Analyzed	Dil Fac
13C3 HFPO-DA	75		50 - 200	02/20/18 10:22	02/21/18 08:24	1

TestAmerica Denver

# Client Sample Results

Client: Chemours Company FC, LLC The  
Project/Site: FAY-OLD OUTFALL SAMP 2/18

TestAmerica Job ID: 280-106447-1

Client Sample ID: FAY-VES-OLDOUTFALLCREEK-A

Lab Sample ID: 280-106447-8

Date Collected: 02/14/18 13:22

Matrix: Water

Date Received: 02/15/18 09:30

## Method: 8321A - HFPO-DA

Analyte	Result	Qualifier	RL	MDL	Unit	D	Prepared	Analyzed	Dil Fac
HFPO-DA	2.3		0.010		ug/L		02/20/18 10:22	02/21/18 09:17	2
Surrogate	%Recovery	Qualifier	Limits				Prepared	Analyzed	Dil Fac
13C3 HFPO-DA	83	D	50 - 200				02/20/18 10:22	02/21/18 09:17	2

TestAmerica Denver

# Client Sample Results

Client: Chemours Company FC, LLC The  
Project/Site: FAY-OLD OUTFALL SAMP 2/18

TestAmerica Job ID: 280-106447-1

Client Sample ID: FAY-VES-OLDOUTFALLCREEK-A2

Lab Sample ID: 280-106447-9

Date Collected: 02/14/18 16:30

Matrix: Water

Date Received: 02/15/18 09:30

## Method: 8321A - HFPO-DA

Analyte	Result	Qualifier	RL	MDL	Unit	D	Prepared	Analyzed	Dil Fac
HFPO-DA	2.5		0.010		ug/L		02/20/18 10:22	02/21/18 09:20	2
Surrogate	%Recovery	Qualifier	Limits				Prepared	Analyzed	Dil Fac
13C3 HFPO-DA	81	D	50 - 200				02/20/18 10:22	02/21/18 09:20	2

TestAmerica Denver

# Client Sample Results

Client: Chemours Company FC, LLC The  
Project/Site: FAY-OLD OUTFALL SAMP 2/18

TestAmerica Job ID: 280-106447-1

Client Sample ID: FAY-VES-OLDOUTFALLCREEK-A3

Lab Sample ID: 280-106447-10

Date Collected: 02/14/18 16:45

Matrix: Water

Date Received: 02/15/18 09:30

## Method: 8321A - HFPO-DA

Analyte	Result	Qualifier	RL	MDL	Unit	D	Prepared	Analyzed	Dil Fac
HFPO-DA	2.2		0.010		ug/L		02/20/18 10:22	02/21/18 09:26	2
Surrogate	%Recovery	Qualifier	Limits				Prepared	Analyzed	Dil Fac
13C3 HFPO-DA	83	D	50 - 200				02/20/18 10:22	02/21/18 09:26	2

TestAmerica Denver



# Client Sample Results

Client: Chemours Company FC, LLC The  
Project/Site: FAY-OLD OUTFALL SAMP 2/18

TestAmerica Job ID: 280-106447-1

Client Sample ID: FAY-VES-OLDOUTFALLCREEKWATERBO

Lab Sample ID: 280-106447-11

Date Collected: 02/14/18 17:00

Matrix: Water

Date Received: 02/15/18 09:30

## Method: 8321A - HFPO-DA

Analyte	Result	Qualifier	RL	MDL	Unit	D	Prepared	Analyzed	Dil Fac
HFPO-DA	1.8		0.010		ug/L		02/20/18 10:22	02/21/18 09:30	2
Surrogate	%Recovery	Qualifier	Limits				Prepared	Analyzed	Dil Fac
13C3 HFPO-DA	85	D	50 - 200				02/20/18 10:22	02/21/18 09:30	2

TestAmerica Denver

# Client Sample Results

Client: Chemours Company FC, LLC The  
Project/Site: FAY-OLD OUTFALL SAMP 2/18

TestAmerica Job ID: 280-106447-1

Client Sample ID: FAY-VES-FB-021418

Lab Sample ID: 280-106447-12

Date Collected: 02/14/18 07:00

Matrix: Water

Date Received: 02/15/18 09:30

## Method: 8321A - HFPO-DA

Analyte	Result	Qualifier	RL	MDL	Unit	D	Prepared	Analyzed	Dil Fac
HFPO-DA	<0.010		0.010		ug/L		02/20/18 10:22	02/21/18 08:41	1

Surrogate	%Recovery	Qualifier	Limits	Prepared	Analyzed	Dil Fac
13C3 HFPO-DA	84		50 - 200	02/20/18 10:22	02/21/18 08:41	1

TestAmerica Denver

## Default Detection Limits

Client: Chemours Company FC, LLC The  
Project/Site: FAY-OLD OUTFALL SAMP 2/18

TestAmerica Job ID: 280-106447-1

Method: 8321A - HFPO-DA

Prep: 3535

Analyte	RL	MDL	Units	Method
HFPO-DA	0.010	0.0051	ug/L	8321A

# Surrogate Summary

Client: Chemours Company FC, LLC The  
Project/Site: FAY-OLD OUTFALL SAMP 2/18

TestAmerica Job ID: 280-106447-1

Method: 8321A - HFPO-DA

Matrix: Water

Prep Type: Total/NA

## Percent Surrogate Recovery (Acceptance Limits)

Lab Sample ID	Client Sample ID	HFPODA (50-200)
280-106447-1	FAY-VES-OLDOUTFALL-A-D	92 D
280-106447-2	FAY-VES-OLDOUTFALL-A	89 D
280-106447-2 DU	FAY-VES-OLDOUTFALL-A	79 D
280-106447-2 MS	FAY-VES-OLDOUTFALL-A	89 D
280-106447-3	FAY-VES-OLDOUTFALL-B	90 D
280-106447-4	FAY-VES-OLDOUTFALL-C	86 D
280-106447-5	FAY-VES-OLDOUTFALL-D	86 D
280-106447-6	FAY-VES-OLDOUTFALL-E	91 D
280-106447-7	FAY-VES-OLDOUTFALLSEEP-	75
280-106447-8	FAY-VES-OLDOUTFALLCREEP	83 D
280-106447-9	FAY-VES-OLDOUTFALLCREEP	81 D
280-106447-10	FAY-VES-OLDOUTFALLCREEP	83 D
280-106447-11	FAY-VES-OLDOUTFALLCREEP	85 D
280-106447-12	FAY-VES-FB-021418	84
DLCK 280-404345/13	Lab Control Sample	104
LCS 280-405473/2-A	Lab Control Sample	81
LCSD 280-405473/4-A	Lab Control Sample Dup	79
LLCS 280-405473/3-A	Lab Control Sample	83
MB 280-405473/1-A	Method Blank	77

### Surrogate Legend

HFPODA = 13C3 HFPO-DA

# QC Sample Results

Client: Chemours Company FC, LLC The  
Project/Site: FAY-OLD OUTFALL SAMP 2/18

TestAmerica Job ID: 280-106447-1

Method: 8321A - HFPO-DA

Lab Sample ID: DLCK 280-404345/13

Matrix: Water

Analysis Batch: 404345

Client Sample ID: Lab Control Sample

Prep Type: Total/NA

Analyte	Spike Added	DLCK Result	DLCK Qualifier	Unit	D	%Rec	% Rec. Limits
HFPO-DA	0.250	<0.50		ug/L		90	70 - 130
Surrogate	DLCK %Recovery	DLCK Qualifier	Limits				
13C3 HFPO-DA	104		50 - 200				

Lab Sample ID: MB 280-405473/1-A

Matrix: Water

Analysis Batch: 405660

Client Sample ID: Method Blank

Prep Type: Total/NA

Prep Batch: 405473

Analyte	MB Result	MB Qualifier	RL	MDL	Unit	D	Prepared	Analyzed	Dil Fac
HFPO-DA	<0.010		0.010		ug/L		02/20/18 10:22	02/21/18 07:42	1
Surrogate	MB %Recovery	MB Qualifier	Limits	Prepared	Analyzed	Dil Fac			
13C3 HFPO-DA	77		50 - 200	02/20/18 10:22	02/21/18 07:42	1			

Lab Sample ID: LCS 280-405473/2-A

Matrix: Water

Analysis Batch: 405660

Client Sample ID: Lab Control Sample

Prep Type: Total/NA

Prep Batch: 405473

Analyte	Spike Added	LCS Result	LCS Qualifier	Unit	D	%Rec	% Rec. Limits
HFPO-DA	0.200	0.191		ug/L		95	70 - 130
Surrogate	LCS %Recovery	LCS Qualifier	Limits				
13C3 HFPO-DA	81		50 - 200				

Lab Sample ID: LCSD 280-405473/4-A

Matrix: Water

Analysis Batch: 405660

Client Sample ID: Lab Control Sample Dup

Prep Type: Total/NA

Prep Batch: 405473

Analyte	Spike Added	LCSD Result	LCSD Qualifier	Unit	D	%Rec	% Rec. Limits	RPD	Limit
HFPO-DA	0.200	0.197		ug/L		99	70 - 130	3	20
Surrogate	LCSD %Recovery	LCSD Qualifier	Limits	RPD	Limit				
13C3 HFPO-DA	79		50 - 200						

Lab Sample ID: LLCS 280-405473/3-A

Matrix: Water

Analysis Batch: 405660

Client Sample ID: Lab Control Sample

Prep Type: Total/NA

Prep Batch: 405473

Analyte	Spike Added	LLCS Result	LLCS Qualifier	Unit	D	%Rec	% Rec. Limits
HFPO-DA	0.0200	0.0193		ug/L		97	70 - 130
Surrogate	LLCS %Recovery	LLCS Qualifier	Limits				
13C3 HFPO-DA	83		50 - 200				

TestAmerica Denver

# QC Sample Results

Client: Chemours Company FC, LLC The  
Project/Site: FAY-OLD OUTFALL SAMP 2/18

TestAmerica Job ID: 280-106447-1

## Method: 8321A - HFPO-DA (Continued)

Lab Sample ID: 280-106447-2 MS

Matrix: Water

Analysis Batch: 405660

Client Sample ID: FAY-VES-OLDOUTFALL-A

Prep Type: Total/NA

Prep Batch: 405473

Analyte	Sample Result	Sample Qualifier	Spike Added	MS Result	MS Qualifier	Unit	D	%Rec	Limits
HFPO-DA	7.4		0.180	7.71	4	ug/L		163	70 - 130
Surrogate	MS %Recovery	MS Qualifier	Limits						
13C3 HFPO-DA	89	D	50 - 200						

Lab Sample ID: 280-106447-2 DU

Matrix: Water

Analysis Batch: 405660

Client Sample ID: FAY-VES-OLDOUTFALL-A

Prep Type: Total/NA

Prep Batch: 405473

Analyte	Sample Result	Sample Qualifier	DU Result	DU Qualifier	Unit	D	RPD	Limit
HFPO-DA	7.4		7.77		ug/L		5	20
Surrogate	DU %Recovery	DU Qualifier	Limits					
13C3 HFPO-DA	79	D	50 - 200					

TestAmerica Denver

# QC Association Summary

Client: Chemours Company FC, LLC The  
Project/Site: FAY-OLD OUTFALL SAMP 2/18

TestAmerica Job ID: 280-106447-1

## LCMS

### Analysis Batch: 404345

Lab Sample ID	Client Sample ID	Prep Type	Matrix	Method	Prep Batch
DLCK 280-404345/13	Lab Control Sample	Total/NA	Water	8321A	

### Prep Batch: 405473

Lab Sample ID	Client Sample ID	Prep Type	Matrix	Method	Prep Batch
280-106447-1	FAY-VES-OLDOUTFALL-A-D	Total/NA	Water	3535	
280-106447-2	FAY-VES-OLDOUTFALL-A	Total/NA	Water	3535	
280-106447-3	FAY-VES-OLDOUTFALL-B	Total/NA	Water	3535	
280-106447-4	FAY-VES-OLDOUTFALL-C	Total/NA	Water	3535	
280-106447-5	FAY-VES-OLDOUTFALL-D	Total/NA	Water	3535	
280-106447-6	FAY-VES-OLDOUTFALL-E	Total/NA	Water	3535	
280-106447-7	FAY-VES-OLDOUTFALLSEEP-A	Total/NA	Water	3535	
280-106447-8	FAY-VES-OLDOUTFALLCREEK-A	Total/NA	Water	3535	
280-106447-9	FAY-VES-OLDOUTFALLCREEK-A2	Total/NA	Water	3535	
280-106447-10	FAY-VES-OLDOUTFALLCREEK-A3	Total/NA	Water	3535	
280-106447-11	FAY-VES-OLDOUTFALLCREEKWATERBO	Total/NA	Water	3535	
280-106447-12	FAY-VES-FB-021418	Total/NA	Water	3535	
MB 280-405473/1-A	Method Blank	Total/NA	Water	3535	
LCS 280-405473/2-A	Lab Control Sample	Total/NA	Water	3535	
LCSD 280-405473/4-A	Lab Control Sample Dup	Total/NA	Water	3535	
LLCS 280-405473/3-A	Lab Control Sample	Total/NA	Water	3535	
280-106447-2 MS	FAY-VES-OLDOUTFALL-A	Total/NA	Water	3535	
280-106447-2 DU	FAY-VES-OLDOUTFALL-A	Total/NA	Water	3535	

### Analysis Batch: 405660

Lab Sample ID	Client Sample ID	Prep Type	Matrix	Method	Prep Batch
280-106447-1	FAY-VES-OLDOUTFALL-A-D	Total/NA	Water	8321A	405473
280-106447-2	FAY-VES-OLDOUTFALL-A	Total/NA	Water	8321A	405473
280-106447-3	FAY-VES-OLDOUTFALL-B	Total/NA	Water	8321A	405473
280-106447-4	FAY-VES-OLDOUTFALL-C	Total/NA	Water	8321A	405473
280-106447-5	FAY-VES-OLDOUTFALL-D	Total/NA	Water	8321A	405473
280-106447-6	FAY-VES-OLDOUTFALL-E	Total/NA	Water	8321A	405473
280-106447-7	FAY-VES-OLDOUTFALLSEEP-A	Total/NA	Water	8321A	405473
280-106447-8	FAY-VES-OLDOUTFALLCREEK-A	Total/NA	Water	8321A	405473
280-106447-9	FAY-VES-OLDOUTFALLCREEK-A2	Total/NA	Water	8321A	405473
280-106447-10	FAY-VES-OLDOUTFALLCREEK-A3	Total/NA	Water	8321A	405473
280-106447-11	FAY-VES-OLDOUTFALLCREEKWATERBO	Total/NA	Water	8321A	405473
280-106447-12	FAY-VES-FB-021418	Total/NA	Water	8321A	405473
MB 280-405473/1-A	Method Blank	Total/NA	Water	8321A	405473
LCS 280-405473/2-A	Lab Control Sample	Total/NA	Water	8321A	405473
LCSD 280-405473/4-A	Lab Control Sample Dup	Total/NA	Water	8321A	405473
LLCS 280-405473/3-A	Lab Control Sample	Total/NA	Water	8321A	405473
280-106447-2 MS	FAY-VES-OLDOUTFALL-A	Total/NA	Water	8321A	405473
280-106447-2 DU	FAY-VES-OLDOUTFALL-A	Total/NA	Water	8321A	405473

# Lab Chronicle

Client: Chemours Company FC, LLC The  
Project/Site: FAY-OLD OUTFALL SAMP 2/18

TestAmerica Job ID: 280-106447-1

Client Sample ID: FAY-VES-OLDOUTFALL-A-D

Lab Sample ID: 280-106447-1

Date Collected: 02/14/18 12:00

Matrix: Water

Date Received: 02/15/18 09:30

Prep Type	Batch Type	Batch Method	Run	Dil Factor	Initial Amount	Final Amount	Batch Number	Prepared or Analyzed	Analyst	Lab
Total/NA	Prep	3535			288.9 mL	5 mL	405473	02/20/18 10:22	DFB1	TAL DEN
Total/NA	Analysis	8321A		10			405660	02/21/18 08:50	AGCM	TAL DEN

Client Sample ID: FAY-VES-OLDOUTFALL-A

Lab Sample ID: 280-106447-2

Date Collected: 02/14/18 12:00

Matrix: Water

Date Received: 02/15/18 09:30

Prep Type	Batch Type	Batch Method	Run	Dil Factor	Initial Amount	Final Amount	Batch Number	Prepared or Analyzed	Analyst	Lab
Total/NA	Prep	3535			295.8 mL	5 mL	405473	02/20/18 10:22	DFB1	TAL DEN
Total/NA	Analysis	8321A		10			405660	02/21/18 08:54	AGCM	TAL DEN

Client Sample ID: FAY-VES-OLDOUTFALL-B

Lab Sample ID: 280-106447-3

Date Collected: 02/14/18 12:45

Matrix: Water

Date Received: 02/15/18 09:30

Prep Type	Batch Type	Batch Method	Run	Dil Factor	Initial Amount	Final Amount	Batch Number	Prepared or Analyzed	Analyst	Lab
Total/NA	Prep	3535			292.6 mL	5 mL	405473	02/20/18 10:22	DFB1	TAL DEN
Total/NA	Analysis	8321A		10			405660	02/21/18 09:04	AGCM	TAL DEN

Client Sample ID: FAY-VES-OLDOUTFALL-C

Lab Sample ID: 280-106447-4

Date Collected: 02/14/18 13:38

Matrix: Water

Date Received: 02/15/18 09:30

Prep Type	Batch Type	Batch Method	Run	Dil Factor	Initial Amount	Final Amount	Batch Number	Prepared or Analyzed	Analyst	Lab
Total/NA	Prep	3535			270.9 mL	5 mL	405473	02/20/18 10:22	DFB1	TAL DEN
Total/NA	Analysis	8321A		10			405660	02/21/18 09:07	AGCM	TAL DEN

Client Sample ID: FAY-VES-OLDOUTFALL-D

Lab Sample ID: 280-106447-5

Date Collected: 02/14/18 14:05

Matrix: Water

Date Received: 02/15/18 09:30

Prep Type	Batch Type	Batch Method	Run	Dil Factor	Initial Amount	Final Amount	Batch Number	Prepared or Analyzed	Analyst	Lab
Total/NA	Prep	3535			283.5 mL	5 mL	405473	02/20/18 10:22	DFB1	TAL DEN
Total/NA	Analysis	8321A		10			405660	02/21/18 09:10	AGCM	TAL DEN

Client Sample ID: FAY-VES-OLDOUTFALL-E

Lab Sample ID: 280-106447-6

Date Collected: 02/14/18 14:15

Matrix: Water

Date Received: 02/15/18 09:30

Prep Type	Batch Type	Batch Method	Run	Dil Factor	Initial Amount	Final Amount	Batch Number	Prepared or Analyzed	Analyst	Lab
Total/NA	Prep	3535			281 mL	5 mL	405473	02/20/18 10:22	DFB1	TAL DEN
Total/NA	Analysis	8321A		10			405660	02/21/18 09:13	AGCM	TAL DEN

TestAmerica Denver



# Lab Chronicle

Client: Chemours Company FC, LLC The  
Project/Site: FAY-OLD OUTFALL SAMP 2/18

TestAmerica Job ID: 280-106447-1

Client Sample ID: FAY-VES-OLDOUTFALLSEEP-A

Lab Sample ID: 280-106447-7

Date Collected: 02/14/18 13:10

Matrix: Water

Date Received: 02/15/18 09:30

Prep Type	Batch Type	Batch Method	Run	Dil Factor	Initial Amount	Final Amount	Batch Number	Prepared or Analyzed	Analyst	Lab
Total/NA	Prep	3535			287.4 mL	5 mL	405473	02/20/18 10:22	DFB1	TAL DEN
Total/NA	Analysis	8321A		1			405660	02/21/18 08:24	AGCM	TAL DEN

Client Sample ID: FAY-VES-OLDOUTFALLCREEK-A

Lab Sample ID: 280-106447-8

Date Collected: 02/14/18 13:22

Matrix: Water

Date Received: 02/15/18 09:30

Prep Type	Batch Type	Batch Method	Run	Dil Factor	Initial Amount	Final Amount	Batch Number	Prepared or Analyzed	Analyst	Lab
Total/NA	Prep	3535			293.6 mL	5 mL	405473	02/20/18 10:22	DFB1	TAL DEN
Total/NA	Analysis	8321A		2			405660	02/21/18 09:17	AGCM	TAL DEN

Client Sample ID: FAY-VES-OLDOUTFALLCREEK-A2

Lab Sample ID: 280-106447-9

Date Collected: 02/14/18 16:30

Matrix: Water

Date Received: 02/15/18 09:30

Prep Type	Batch Type	Batch Method	Run	Dil Factor	Initial Amount	Final Amount	Batch Number	Prepared or Analyzed	Analyst	Lab
Total/NA	Prep	3535			290 mL	5 mL	405473	02/20/18 10:22	DFB1	TAL DEN
Total/NA	Analysis	8321A		2			405660	02/21/18 09:20	AGCM	TAL DEN

Client Sample ID: FAY-VES-OLDOUTFALLCREEK-A3

Lab Sample ID: 280-106447-10

Date Collected: 02/14/18 16:45

Matrix: Water

Date Received: 02/15/18 09:30

Prep Type	Batch Type	Batch Method	Run	Dil Factor	Initial Amount	Final Amount	Batch Number	Prepared or Analyzed	Analyst	Lab
Total/NA	Prep	3535			294.9 mL	5 mL	405473	02/20/18 10:22	DFB1	TAL DEN
Total/NA	Analysis	8321A		2			405660	02/21/18 09:26	AGCM	TAL DEN

Client Sample ID: FAY-VES-OLDOUTFALLCREEKWATERBO

Lab Sample ID: 280-106447-11

Date Collected: 02/14/18 17:00

Matrix: Water

Date Received: 02/15/18 09:30

Prep Type	Batch Type	Batch Method	Run	Dil Factor	Initial Amount	Final Amount	Batch Number	Prepared or Analyzed	Analyst	Lab
Total/NA	Prep	3535			294.1 mL	5 mL	405473	02/20/18 10:22	DFB1	TAL DEN
Total/NA	Analysis	8321A		2			405660	02/21/18 09:30	AGCM	TAL DEN

Client Sample ID: FAY-VES-FB-021418

Lab Sample ID: 280-106447-12

Date Collected: 02/14/18 07:00

Matrix: Water

Date Received: 02/15/18 09:30

Prep Type	Batch Type	Batch Method	Run	Dil Factor	Initial Amount	Final Amount	Batch Number	Prepared or Analyzed	Analyst	Lab
Total/NA	Prep	3535			293.5 mL	5 mL	405473	02/20/18 10:22	DFB1	TAL DEN
Total/NA	Analysis	8321A		1			405660	02/21/18 08:41	AGCM	TAL DEN

TestAmerica Denver

# Lab Chronicle

Client: Chemours Company FC, LLC The  
Project/Site: FAY-OLD OUTFALL SAMP 2/18

TestAmerica Job ID: 280-106447-1

Client Sample ID: Method Blank

Lab Sample ID: MB 280-405473/1-A

Date Collected: N/A

Matrix: Water

Date Received: N/A

Prep Type	Batch Type	Batch Method	Run	Dil Factor	Initial Amount	Final Amount	Batch Number	Prepared or Analyzed	Analyst	Lab
Total/NA	Prep	3535			250 mL	5 mL	405473	02/20/18 10:22	DFB1	TAL DEN
Total/NA	Analysis	8321A		1			405660	02/21/18 07:42	AGCM	TAL DEN

Client Sample ID: Lab Control Sample

Lab Sample ID: DLCK 280-404345/13

Date Collected: N/A

Matrix: Water

Date Received: N/A

Prep Type	Batch Type	Batch Method	Run	Dil Factor	Initial Amount	Final Amount	Batch Number	Prepared or Analyzed	Analyst	Lab
Total/NA	Analysis	8321A		1			404345	02/08/18 13:38	AGCM	TAL DEN

Client Sample ID: Lab Control Sample

Lab Sample ID: LCS 280-405473/2-A

Date Collected: N/A

Matrix: Water

Date Received: N/A

Prep Type	Batch Type	Batch Method	Run	Dil Factor	Initial Amount	Final Amount	Batch Number	Prepared or Analyzed	Analyst	Lab
Total/NA	Prep	3535			250 mL	5 mL	405473	02/20/18 10:22	DFB1	TAL DEN
Total/NA	Analysis	8321A		1			405660	02/21/18 07:45	AGCM	TAL DEN

Client Sample ID: Lab Control Sample Dup

Lab Sample ID: LCSD 280-405473/4-A

Date Collected: N/A

Matrix: Water

Date Received: N/A

Prep Type	Batch Type	Batch Method	Run	Dil Factor	Initial Amount	Final Amount	Batch Number	Prepared or Analyzed	Analyst	Lab
Total/NA	Prep	3535			250 mL	5 mL	405473	02/20/18 10:22	DFB1	TAL DEN
Total/NA	Analysis	8321A		1			405660	02/21/18 07:52	AGCM	TAL DEN

Client Sample ID: Lab Control Sample

Lab Sample ID: LLCS 280-405473/3-A

Date Collected: N/A

Matrix: Water

Date Received: N/A

Prep Type	Batch Type	Batch Method	Run	Dil Factor	Initial Amount	Final Amount	Batch Number	Prepared or Analyzed	Analyst	Lab
Total/NA	Prep	3535			250 mL	5 mL	405473	02/20/18 10:22	DFB1	TAL DEN
Total/NA	Analysis	8321A		1			405660	02/21/18 07:49	AGCM	TAL DEN

Client Sample ID: FAY-VES-OLDOUTFALL-A

Lab Sample ID: 280-106447-2 MS

Date Collected: 02/14/18 12:00

Matrix: Water

Date Received: 02/15/18 09:30

Prep Type	Batch Type	Batch Method	Run	Dil Factor	Initial Amount	Final Amount	Batch Number	Prepared or Analyzed	Analyst	Lab
Total/NA	Prep	3535			278.2 mL	5 mL	405473	02/20/18 10:22	DFB1	TAL DEN
Total/NA	Analysis	8321A		10			405660	02/21/18 09:00	AGCM	TAL DEN

TestAmerica Denver

# Lab Chronicle

Client: Chemours Company FC, LLC The  
Project/Site: FAY-OLD OUTFALL SAMP 2/18

TestAmerica Job ID: 280-106447-1

Client Sample ID: FAY-VES-OLDOUTFALL-A

Lab Sample ID: 280-106447-2 DU

Date Collected: 02/14/18 12:00

Matrix: Water

Date Received: 02/15/18 09:30

Prep Type	Batch Type	Batch Method	Run	Dil Factor	Initial Amount	Final Amount	Batch Number	Prepared or Analyzed	Analyst	Lab
Total/NA	Prep	3535			280.9 mL	5 mL	405473	02/20/18 10:22	DFB1	TAL DEN
Total/NA	Analysis	8321A		10			405660	02/21/18 08:57	AGCM	TAL DEN

## Laboratory References:

TAL DEN = TestAmerica Denver, 4955 Yarrow Street, Arvada, CO 80002, TEL (303)736-0100

## Accreditation/Certification Summary

Client: Chemours Company FC, LLC The  
Project/Site: FAY-OLD OUTFALL SAMP 2/18

TestAmerica Job ID: 280-106447-1

### Laboratory: TestAmerica Denver

Unless otherwise noted, all analytes for this laboratory were covered under each accreditation/certification below.

Authority	Program	EPA Region	Identification Number	Expiration Date
North Carolina (WW/SW)	State Program	4	358	12-31-18

The following analytes are included in this report, but accreditation/certification is not offered by the governing authority:

Analysis Method	Prep Method	Matrix	Analyte
8321A	3535	Water	HFPO-DA

## Method Summary

Client: Chemours Company FC, LLC The  
Project/Site: FAY-OLD OUTFALL SAMP 2/18

TestAmerica Job ID: 280-106447-1

---

---

Method	Method Description	Protocol	Laboratory
8321A	HFPO-DA	SW846	TAL DEN

---

### Protocol References:

SW846 = "Test Methods For Evaluating Solid Waste, Physical/Chemical Methods", Third Edition, November 1986 And Its Updates.

### Laboratory References:

TAL DEN = TestAmerica Denver, 4955 Yarrow Street, Arvada, CO 80002, TEL (303)736-0100

## Sample Summary

Client: Chemours Company FC, LLC The  
Project/Site: FAY-OLD OUTFALL SAMP 2/18

TestAmerica Job ID: 280-106447-1

Lab Sample ID	Client Sample ID	Matrix	Collected	Received
280-106447-1	FAY-VES-OLDOUTFALL-A-D	Water	02/14/18 12:00	02/15/18 09:30
280-106447-2	FAY-VES-OLDOUTFALL-A	Water	02/14/18 12:00	02/15/18 09:30
280-106447-3	FAY-VES-OLDOUTFALL-B	Water	02/14/18 12:45	02/15/18 09:30
280-106447-4	FAY-VES-OLDOUTFALL-C	Water	02/14/18 13:38	02/15/18 09:30
280-106447-5	FAY-VES-OLDOUTFALL-D	Water	02/14/18 14:05	02/15/18 09:30
280-106447-6	FAY-VES-OLDOUTFALL-E	Water	02/14/18 14:15	02/15/18 09:30
280-106447-7	FAY-VES-OLDOUTFALLSEEP-A	Water	02/14/18 13:10	02/15/18 09:30
280-106447-8	FAY-VES-OLDOUTFALLCREEK-A	Water	02/14/18 13:22	02/15/18 09:30
280-106447-9	FAY-VES-OLDOUTFALLCREEK-A2	Water	02/14/18 16:30	02/15/18 09:30
280-106447-10	FAY-VES-OLDOUTFALLCREEK-A3	Water	02/14/18 16:45	02/15/18 09:30
280-106447-11	FAY-VES-OLDOUTFALLCREEKWATERBO	Water	02/14/18 17:00	02/15/18 09:30
280-106447-12	FAY-VES-FB-021418	Water	02/14/18 07:00	02/15/18 09:30

TestAmerica Denver

## LCMS MANUAL INTEGRATION SUMMARY

Lab Name: TestAmerica Denver Job No.: 280-106447-1

SDG No.: \_\_\_\_\_

Instrument ID: LC\_LCMS7 Analysis Batch Number: 404345Lab Sample ID: STD001 280-404345/3 IC Client Sample ID: \_\_\_\_\_Date Analyzed: 02/08/18 13:05 Lab File ID: hfpo718B08034.d GC Column: Synergi Hydro ID: \_\_\_\_\_

COMPOUND NAME	RETENTION TIME	MANUAL INTEGRATION		
		REASON	ANALYST	DATE
HFPO-DA	1.06	Assign Peak	meyera	02/08/18 15:19

Lab Sample ID: STD002 280-404345/4 IC Client Sample ID: \_\_\_\_\_Date Analyzed: 02/08/18 13:08 Lab File ID: hfpo718B08035.d GC Column: Synergi Hydro ID: \_\_\_\_\_

COMPOUND NAME	RETENTION TIME	MANUAL INTEGRATION		
		REASON	ANALYST	DATE
HFPO-DA	1.06	Baseline	meyera	02/08/18 15:19

Lab Sample ID: DLCK 280-404345/13 Client Sample ID: \_\_\_\_\_Date Analyzed: 02/08/18 13:38 Lab File ID: hfpo718B08044.d GC Column: Synergi Hydro ID: \_\_\_\_\_

COMPOUND NAME	RETENTION TIME	MANUAL INTEGRATION		
		REASON	ANALYST	DATE
HFPO-DA	1.06	Baseline	meyera	02/08/18 15:20

# REAGENT TRACEABILITY SUMMARY

Lab Name: TestAmerica Denver

Job No.: 280-106447-1

SDG No.: \_\_\_\_\_

Reagent ID	Exp Date	Prep Date	Dilutant Used	Reagent Final Volume	Parent Reagent		Analyte	Concentration
					Reagent ID	Volume Added		
<b>HFPO I.S._00008</b>	12/12/18	01/30/18	LCMS Grade MeOH, Lot LCMS_MeOH_00110	100 mL	13C3 HFPO-DA_00008	1 mL	13C3 HFPO-DA	0.5 ug/mL
							13C3 HFPO-DA (IS)	0.5 ug/mL
.13C3 HFPO-DA_00008	01/30/19	Wellington Laboratories, Lot M3HFPOADA0817			(Purchased Reagent)		13C3 HFPO-DA	50 ug/mL
							13C3 HFPO-DA (IS)	50 ug/mL
<b>HFPO Spike_00004</b>	10/30/18	10/30/17	LCMS Grade MeOH, Lot LCMS MeOH 00110	100 mL	HFPO-DA_00004	1 mL	HFPO-DA	0.5 ug/mL
.HFPO-DA_00004	07/13/20	Wellington Laboratories, Lot HFPODA0717			(Purchased Reagent)		HFPO-DA	50 ug/mL
<b>HFPO_CAL-0_00032</b>	02/22/18	02/08/18	PFC Dill_Solvent, Lot 00016	1 mL	HFPO I.S._00008	20 uL	13C3 HFPO-DA	10 ug/L
.HFPO I.S._00008	12/12/18	01/30/18	LCMS Grade MeOH, Lot LCMS MeOH 00110	100 mL	13C3 HFPO-DA_00008	1 mL	13C3 HFPO-DA	0.5 ug/mL
..13C3 HFPO-DA_00008	01/30/19	Wellington Laboratories, Lot M3HFPOADA0817			(Purchased Reagent)		13C3 HFPO-DA	50 ug/mL
<b>HFPO_CAL-1_00032</b>	02/22/18	02/08/18	80:20 Methanol : H2O, Lot 00016	1 mL	HFPO I.S._00008	20 uL	13C3 HFPO-DA	10 ug/L
							13C3 HFPO-DA (IS)	10 ug/L
					HFPO Spike_00004	0.5 uL	HFPO-DA	0.25 ug/L
.HFPO I.S._00008	12/12/18	01/30/18	LCMS Grade MeOH, Lot LCMS_MeOH_00110	100 mL	13C3 HFPO-DA_00008	1 mL	13C3 HFPO-DA	0.5 ug/mL
							13C3 HFPO-DA (IS)	0.5 ug/mL
..13C3 HFPO-DA_00008	01/30/19	Wellington Laboratories, Lot M3HFPOADA0817			(Purchased Reagent)		13C3 HFPO-DA	50 ug/mL
							13C3 HFPO-DA (IS)	50 ug/mL
.HFPO Spike_00004	10/30/18	10/30/17	LCMS Grade MeOH, Lot LCMS MeOH 00110	100 mL	HFPO-DA_00004	1 mL	HFPO-DA	0.5 ug/mL
..HFPO-DA_00004	07/13/20	Wellington Laboratories, Lot HFPODA0717			(Purchased Reagent)		HFPO-DA	50 ug/mL
<b>HFPO_CAL-2_00033</b>	02/22/18	02/08/18	80:20 Methanol : H2O, Lot 00016	1 mL	HFPO I.S._00008	20 uL	13C3 HFPO-DA	10 ug/L
							13C3 HFPO-DA (IS)	10 ug/L
					HFPO Spike_00004	1 uL	HFPO-DA	0.5 ug/L
.HFPO I.S._00008	12/12/18	01/30/18	LCMS Grade MeOH, Lot LCMS_MeOH_00110	100 mL	13C3 HFPO-DA_00008	1 mL	13C3 HFPO-DA	0.5 ug/mL
							13C3 HFPO-DA (IS)	0.5 ug/mL
..13C3 HFPO-DA_00008	01/30/19	Wellington Laboratories, Lot M3HFPOADA0817			(Purchased Reagent)		13C3 HFPO-DA	50 ug/mL
							13C3 HFPO-DA (IS)	50 ug/mL
.HFPO Spike_00004	10/30/18	10/30/17	LCMS Grade MeOH, Lot LCMS MeOH 00110	100 mL	HFPO-DA_00004	1 mL	HFPO-DA	0.5 ug/mL
..HFPO-DA_00004	07/13/20	Wellington Laboratories, Lot HFPODA0717			(Purchased Reagent)		HFPO-DA	50 ug/mL
<b>HFPO_CAL-3_00032</b>	02/22/18	02/08/18	80:20 Methanol : H2O, Lot 00016	1 mL	HFPO I.S._00008	20 uL	13C3 HFPO-DA	10 ug/L
							13C3 HFPO-DA (IS)	10 ug/L
					HFPO Spike_00004	2 uL	HFPO-DA	1 ug/L
.HFPO I.S._00008	12/12/18	01/30/18	LCMS Grade MeOH, Lot LCMS_MeOH_00110	100 mL	13C3 HFPO-DA_00008	1 mL	13C3 HFPO-DA	0.5 ug/mL
							13C3 HFPO-DA (IS)	0.5 ug/mL



# REAGENT TRACEABILITY SUMMARY

Lab Name: TestAmerica Denver

Job No.: 280-106447-1

SDG No.: \_\_\_\_\_

Reagent ID	Exp Date	Prep Date	Dilutant Used	Reagent Final Volume	Parent Reagent		Analyte	Concentration
					Reagent ID	Volume Added		
..13C3 HFPO-DA_00008	01/30/19		Wellington Laboratories, Lot M3HFPOADA0817		(Purchased Reagent)		13C3 HFPO-DA	50 ug/mL
.HFPO Spike_00004	10/30/18	10/30/17	LCMS Grade MeOH, Lot LCMS MeOH 00110	100 mL	HFPO-DA_00004	1 mL	13C3 HFPO-DA (IS) HFPO-DA	50 ug/mL 0.5 ug/mL
..HFPO-DA_00004	07/13/20		Wellington Laboratories, Lot HFPODA0717		(Purchased Reagent)		HFPO-DA	50 ug/mL
HFPO_CAL-4_00032	02/22/18	02/08/18	80:20 Methanol : H2O, Lot 00016	1 mL	HFPO I.S._00008	20 uL	13C3 HFPO-DA	10 ug/L
					HFPO Spike_00004	4 uL	13C3 HFPO-DA (IS) HFPO-DA	10 ug/L 2 ug/L
.HFPO I.S._00008	12/12/18	01/30/18	LCMS Grade MeOH, Lot LCMS MeOH 00110	100 mL	13C3 HFPO-DA_00008	1 mL	13C3 HFPO-DA	0.5 ug/mL
..13C3 HFPO-DA_00008	01/30/19		Wellington Laboratories, Lot M3HFPOADA0817		(Purchased Reagent)		13C3 HFPO-DA (IS) 13C3 HFPO-DA	0.5 ug/mL 50 ug/mL
.HFPO Spike_00004	10/30/18	10/30/17	LCMS Grade MeOH, Lot LCMS MeOH 00110	100 mL	HFPO-DA_00004	1 mL	13C3 HFPO-DA (IS) HFPO-DA	50 ug/mL 0.5 ug/mL
..HFPO-DA_00004	07/13/20		Wellington Laboratories, Lot HFPODA0717		(Purchased Reagent)		HFPO-DA	50 ug/mL
HFPO_CAL-5_00080	02/22/18	02/08/18	80:20 Methanol : H2O, Lot 00016	1 mL	HFPO I.S._00008	20 uL	13C3 HFPO-DA	10 ug/L
					HFPO Spike_00004	10 uL	13C3 HFPO-DA (IS) HFPO-DA	10 ug/L 5 ug/L
.HFPO I.S._00008	12/12/18	01/30/18	LCMS Grade MeOH, Lot LCMS MeOH 00110	100 mL	13C3 HFPO-DA_00008	1 mL	13C3 HFPO-DA	0.5 ug/mL
..13C3 HFPO-DA_00008	01/30/19		Wellington Laboratories, Lot M3HFPOADA0817		(Purchased Reagent)		13C3 HFPO-DA (IS) 13C3 HFPO-DA	0.5 ug/mL 50 ug/mL
.HFPO Spike_00004	10/30/18	10/30/17	LCMS Grade MeOH, Lot LCMS MeOH 00110	100 mL	HFPO-DA_00004	1 mL	13C3 HFPO-DA (IS) HFPO-DA	50 ug/mL 0.5 ug/mL
..HFPO-DA_00004	07/13/20		Wellington Laboratories, Lot HFPODA0717		(Purchased Reagent)		HFPO-DA	50 ug/mL
HFPO_CAL-5_00081	02/22/18	02/20/18	80:20 Methanol : H2O, Lot 00016	1 mL	HFPO I.S._00009	20 uL	13C3 HFPO-DA	10 ug/L
					HFPO Spike_00004	10 uL	HFPO-DA	5 ug/L
.HFPO I.S._00009	02/20/19	02/20/18	LCMS Grade MeOH, Lot LCMS MeOH 00110	100 mL	13C3 HFPO-DA_00009	1 mL	13C3 HFPO-DA	0.5 ug/mL
..13C3 HFPO-DA_00009	02/20/19		Wellington Laboratories, Lot M3HFPOADA0817		(Purchased Reagent)		13C3 HFPO-DA	50 ug/mL
.HFPO Spike_00004	10/30/18	10/30/17	LCMS Grade MeOH, Lot LCMS MeOH 00110	100 mL	HFPO-DA_00004	1 mL	HFPO-DA	0.5 ug/mL
..HFPO-DA_00004	07/13/20		Wellington Laboratories, Lot HFPODA0717		(Purchased Reagent)		HFPO-DA	50 ug/mL
HFPO_CAL-6_00080	02/22/18	02/08/18	80:20 Methanol : H2O, Lot 00016	1 mL	HFPO I.S._00008	20 uL	13C3 HFPO-DA	10 ug/L
					HFPO Spike_00004	20 uL	13C3 HFPO-DA (IS) HFPO-DA	10 ug/L 10 ug/L
.HFPO I.S._00008	12/12/18	01/30/18	LCMS Grade MeOH, Lot LCMS MeOH 00110	100 mL	13C3 HFPO-DA_00008	1 mL	13C3 HFPO-DA	0.5 ug/mL

REAGENT TRACEABILITY SUMMARY

Lab Name: TestAmerica Denver

Job No.: 280-106447-1

SDG No.: \_\_\_\_\_

Reagent ID	Exp Date	Prep Date	Dilutant Used	Reagent Final Volume	Parent Reagent		Analyte	Concentration
					Reagent ID	Volume Added		
..13C3 HFPO-DA_00008	01/30/19		Wellington Laboratories, Lot M3HFPOADA0817		(Purchased Reagent)		13C3 HFPO-DA (IS)	0.5 ug/mL
							13C3 HFPO-DA	50 ug/mL
.HFPO Spike_00004	10/30/18	10/30/17	LCMS Grade MeOH, Lot LCMS MeOH 00110	100 mL	HFPO-DA_00004	1 mL	13C3 HFPO-DA (IS)	50 ug/mL
							HFPO-DA	0.5 ug/mL
..HFPO-DA_00004	07/13/20		Wellington Laboratories, Lot HFPODA0717		(Purchased Reagent)		HFPO-DA	50 ug/mL
<b>HFPO_CAL-6_00081</b>	02/22/18	02/20/18	80:20 Methanol : H2O, Lot 00016	1 mL	HFPO I.S._00009	20 uL	13C3 HFPO-DA	10 ug/L
					HFPO Spike_00004	20 uL	HFPO-DA	10 ug/L
.HFPO I.S._00009	02/20/19	02/20/18	LCMS Grade MeOH, Lot LCMS MeOH 00110	100 mL	13C3 HFPO-DA_00009	1 mL	13C3 HFPO-DA	0.5 ug/mL
..13C3 HFPO-DA_00009	02/20/19		Wellington Laboratories, Lot M3HFPOADA0817		(Purchased Reagent)		13C3 HFPO-DA	50 ug/mL
.HFPO Spike_00004	10/30/18	10/30/17	LCMS Grade MeOH, Lot LCMS MeOH 00110	100 mL	HFPO-DA_00004	1 mL	HFPO-DA	0.5 ug/mL
..HFPO-DA_00004	07/13/20		Wellington Laboratories, Lot HFPODA0717		(Purchased Reagent)		HFPO-DA	50 ug/mL
<b>HFPO_CAL-7_00032</b>	02/22/18	02/08/18	80:20 Methanol : H2O, Lot 00016	1 mL	HFPO I.S._00008	20 uL	13C3 HFPO-DA	10 ug/L
							13C3 HFPO-DA (IS)	10 ug/L
					HFPO Spike_00004	50 uL	HFPO-DA	25 ug/L
.HFPO I.S._00008	12/12/18	01/30/18	LCMS Grade MeOH, Lot LCMS MeOH 00110	100 mL	13C3 HFPO-DA_00008	1 mL	13C3 HFPO-DA	0.5 ug/mL
..13C3 HFPO-DA_00008	01/30/19		Wellington Laboratories, Lot M3HFPOADA0817		(Purchased Reagent)		13C3 HFPO-DA (IS)	0.5 ug/mL
							13C3 HFPO-DA	50 ug/mL
.HFPO Spike_00004	10/30/18	10/30/17	LCMS Grade MeOH, Lot LCMS MeOH 00110	100 mL	HFPO-DA_00004	1 mL	HFPO-DA	0.5 ug/mL
..HFPO-DA_00004	07/13/20		Wellington Laboratories, Lot HFPODA0717		(Purchased Reagent)		HFPO-DA	50 ug/mL
<b>HFPO_CAL-8_00032</b>	02/22/18	02/08/18	80:20 Methanol : H2O, Lot 00016	1 mL	HFPO I.S._00008	20 uL	13C3 HFPO-DA	10 ug/L
							13C3 HFPO-DA (IS)	10 ug/L
					HFPO Spike_00004	100 uL	HFPO-DA	50 ug/L
.HFPO I.S._00008	12/12/18	01/30/18	LCMS Grade MeOH, Lot LCMS MeOH 00110	100 mL	13C3 HFPO-DA_00008	1 mL	13C3 HFPO-DA	0.5 ug/mL
..13C3 HFPO-DA_00008	01/30/19		Wellington Laboratories, Lot M3HFPOADA0817		(Purchased Reagent)		13C3 HFPO-DA (IS)	0.5 ug/mL
							13C3 HFPO-DA	50 ug/mL
.HFPO Spike_00004	10/30/18	10/30/17	LCMS Grade MeOH, Lot LCMS MeOH 00110	100 mL	HFPO-DA_00004	1 mL	13C3 HFPO-DA (IS)	50 ug/mL
..HFPO-DA_00004	07/13/20		Wellington Laboratories, Lot HFPODA0717		(Purchased Reagent)		HFPO-DA	0.5 ug/mL
							HFPO-DA	50 ug/mL
<b>HFPO_CAL-9_00001</b>	02/22/18	02/08/18	80:20 Methanol : H2O, Lot 00016	1 mL	HFPO I.S._00008	20 uL	13C3 HFPO-DA	10 ug/L
							13C3 HFPO-DA (IS)	10 ug/L
					HFPO Spike_00004	200 uL	HFPO-DA	100 ug/L

# REAGENT TRACEABILITY SUMMARY

Lab Name: TestAmerica Denver

Job No.: 280-106447-1

SDG No.: \_\_\_\_\_

Reagent ID	Exp Date	Prep Date	Dilutant Used	Reagent Final Volume	Parent Reagent		Analyte	Concentration
					Reagent ID	Volume Added		
.HFPO I.S._00008	12/12/18	01/30/18	LCMS Grade MeOH, Lot LCMS_MeOH_00110	100 mL	13C3 HFPO-DA_00008	1 mL	13C3 HFPO-DA	0.5 ug/mL
..13C3 HFPO-DA_00008	01/30/19	Wellington Laboratories, Lot M3HFPOADA0817			(Purchased Reagent)		13C3 HFPO-DA (IS)	0.5 ug/mL
							13C3 HFPO-DA	50 ug/mL
							13C3 HFPO-DA (IS)	50 ug/mL
.HFPO Spike_00004	10/30/18	10/30/17	LCMS Grade MeOH, Lot LCMS MeOH 00110	100 mL	HFPO-DA_00004	1 mL	HFPO-DA	0.5 ug/mL
..HFPO-DA_00004	07/13/20	Wellington Laboratories, Lot HFPODA0717			(Purchased Reagent)		HFPO-DA	50 ug/mL
HFPO_ICV_00034	02/22/18	02/08/18	80:20 Methanol : H2O, Lot 00016	1 mL	HFPO I.S._00008	20 uL	13C3 HFPO-DA	10 ug/L
					HFPO ICV_00001	10 uL	HFPO-DA	1.95009 ug/L
.HFPO I.S._00008	12/12/18	01/30/18	LCMS Grade MeOH, Lot LCMS MeOH 00110	100 mL	13C3 HFPO-DA_00008	1 mL	13C3 HFPO-DA	0.5 ug/mL
..13C3 HFPO-DA_00008	01/30/19	Wellington Laboratories, Lot M3HFPOADA0817			(Purchased Reagent)		13C3 HFPO-DA	50 ug/mL
.HFPO ICV_00001	11/03/18	11/03/17	Methanol, Lot 12345	100 mL	HFPO SS stock_00002	20 uL	HFPO-DA	0.195009 ug/mL
..HFPO SS stock_00002	11/03/18	11/03/17	Methanol, Lot 12345	500 mL	HFPO SS_00003	0.5026 g	HFPO-DA	975.044 ug/mL
...HFPO SS_00003	05/23/21	Synquest Laboratories, Lot Q141-128			(Purchased Reagent)		HFPO-DA	97 %

Reagent

---

**13C3 HFPO-DA\_00008**



# WELLINGTON LABORATORIES

## CERTIFICATE OF ANALYSIS DOCUMENTATION

**PRODUCT CODE:**


M3HFPO-DA

**LOT NUMBER:**

M3HFPODA0817

**COMPOUND:**2,3,3,3-Tetrafluoro-2-(1,1,2,2,3,3,3-heptafluoropropoxy)-<sup>13</sup>C<sub>3</sub>-propanoic acid**STRUCTURE:****CAS #:**

Not available

**MOLECULAR FORMULA:** $^{13}\text{C}_3^{18}\text{O}_3\text{HF}_{11}\text{O}_3$ **CONCENTRATION:**

50 ± 2.5 µg/ml

**CHEMICAL PURITY:**

&gt;98%

**LAST TESTED:** (mm/dd/yyyy)

08/17/2017

**EXPIRY DATE:** (mm/dd/yyyy)

08/17/2020

**RECOMMENDED STORAGE:**

Store ampoule in a cool, dark place

**MOLECULAR WEIGHT:**

333.03

**SOLVENT(S):**

Methanol

**ISOTOPIC PURITY:**≥99% <sup>13</sup>C(<sup>13</sup>C<sub>3</sub>)**DOCUMENTATION/ DATA ATTACHED:**

Figure 1: LC/MS Data (TIC and Mass Spectrum)

Figure 2: LC/MS/MS Data (Selected MRM Transitions)

**ADDITIONAL INFORMATION:**

- See page 2 for further details.
- Contains ~ 1.5% of two constitutional isomers.
- Product is commercially known as GenX.

**FOR LABORATORY USE ONLY: NOT FOR HUMAN OR DRUG USE**

Certified By:

  
B.G. Chittim, General Manager

Date: 08/25/2017

(mm/dd/yyyy)

Wellington Laboratories Inc., 345 Southgate Dr. Guelph ON N1G 3M5 CANADA  
519-822-2436 • Fax: 519-822-2849 • Info@well-labs.com

### INTENDED USE:

The products prepared by Wellington Laboratories Inc. are for laboratory use only. This certified reference material (CRM) was designed to be used as a standard for the identification and/or quantification of the specific chemical compound it contains.

### HAZARDS:

This product should only be used by qualified personnel familiar with its potential hazards and trained in the handling of hazardous chemicals. Due care should be exercised to prevent unnecessary human contact or ingestion. All procedures should be carried out in a well-functioning fume hood and suitable gloves, eye protection, and clothing should be worn at all times. Waste should be disposed of according to national and regional regulations. Safety Data Sheets (SDSs) are available upon request.

### SYNTHESIS / CHARACTERIZATION:

Where possible, all of our products are synthesized using single-product unambiguous routes. They are then characterized, and their structures and purities confirmed, using a combination of the most relevant techniques, such as NMR, GC/MS, LC/MS/MS, SFC/UV/MS/MS, x-ray crystallography, and melting point. Isotopic purities of mass-labelled compounds are also confirmed using HRGC/HRMS and/or LC/MS/MS.

### HOMOGENEITY:

Prior to solution preparation, crystalline material is tested for homogeneity using a variety of techniques (as stated above) and its solubility in a given diluent is taken into consideration. Duplicate solutions of a new product are prepared from the same crystalline lot and, after the addition of an appropriate internal standard, they are compared by GC/MS, LC/MS/MS and/or SFC/UV/MS/MS. The relative response factors of the analyte of interest in each solution are required to be <5% RSD. New solution lots of existing products are compared to older lots in the same manner, which further confirms the homogeneity of the crystalline material as well as the stability and homogeneity of the solutions in the storage containers. In order to maintain the integrity of the assigned value(s), and associated uncertainty, the dilution or injection of a subsample of this product should be performed using calibrated measuring equipment.

### UNCERTAINTY:

The maximum combined relative standard uncertainty of our reference standard solutions is calculated using the following equation:

The combined relative standard uncertainty,  $u_c(y)$ , of a value  $y$  and the uncertainty of the independent parameters

$x_1, x_2, \dots, x_n$  on which it depends is:

$$u_c(y(x_1, x_2, \dots, x_n)) = \sqrt{\sum_{i=1}^n u(y, x_i)^2}$$

where  $x$  is expressed as a relative standard uncertainty of the individual parameter.

The individual uncertainties taken into account include those associated with weights (calibration of the balance) and volumes (calibration of the volumetric glassware). An expanded maximum combined percent relative uncertainty of  $\pm 5\%$  (calculated with a coverage factor of 2 and a level of confidence of 95%) is stated on the Certificate of Analysis for all of our products.

### TRACEABILITY:

All reference standard solutions are traceable to specific crystalline lots. The microbalances used for solution preparation are regularly tested by an external ISO/IEC 17025 accredited calibration company. In addition, their calibration is verified prior to each weighing using calibrated NIST and/or NRC traceable external weights. All volumetric glassware used is calibrated, of Class A tolerance, and has been tested according to the appropriate ASTM procedures, which are ultimately traceable to NIST. For certain products, traceability to international interlaboratory studies has also been established.

### EXPIRY DATE / PERIOD OF VALIDITY:

Ongoing stability studies of this product have demonstrated stability in its composition and concentration, until the specified expiry date, in the unopened ampoule. Monitoring for any degradation or change in concentration of the listed analyte(s) is performed on a routine basis.

### LIMITED WARRANTY:

At the time of shipment, all products are warranted to be free of defects in material and workmanship and to conform to the stated technical and purity specifications.



### QUALITY MANAGEMENT:

This product was produced using a Quality Management System registered to the latest versions of ISO 9001 by SAI Global, ISO/IEC 17025 by the Canadian Association for Laboratory Accreditation Inc. (CALA; A 1226), and ISO GUIDE 34 by ANSI-ASQ National Accreditation Board (ANAB; AR-1523).



\*\*For additional information or assistance concerning this or any other products from Wellington Laboratories Inc., please visit our website at [www.well-labs.com](http://www.well-labs.com) or contact us directly at [info@well-labs.com](mailto:info@well-labs.com)\*\*

**Figure 1: M3HFPO-DA; LC/MS Data (TIC and Mass Spectrum)**



**Conditions for Figure 1:**

**LC:** Waters Acquity Ultra Performance LC  
**MS:** Micromass Quattro *micro* API MS

**Chromatographic Conditions**

**Column:** Acquity UPLC BEH Shield RP<sub>18</sub>  
1.7  $\mu$ m, 2.1 x 100 mm

**Mobile phase:** Gradient  
Start: 55% MeOH / 45% H<sub>2</sub>O with 10 mM NH<sub>4</sub>OAc buffer  
Ramp to 90% organic over 7.5 min and hold for 1.5 min  
before returning to initial conditions in 0.5 min.

**Time:** 10 min


**Flow:** 300  $\mu$ l/min

**MS Parameters**

**Experiment:** Full Scan (225 - 850 amu)

**Source:** Electrospray (negative)  
**Capillary Voltage (kV)** = 3.00  
**Cone Voltage (V)** = 10.00  
**Cone Gas Flow (l/hr)** = 100  
**Desolvation Gas Flow (l/hr)** = 750

**Figure 2: M3HFPO-DA; LC/MS/MS Data (Selected MRM Transitions)**



**Conditions for Figure 2:**

Injection: Direct loop Injection  
10  $\mu\text{l}$  (500 ng/ml M3HFPO-DA)

Mobile phase: Isocratic 80% MeOH / 20%  $\text{H}_2\text{O}$  with 10 mM  $\text{NH}_4\text{OAc}$  buffer

Flow: 300  $\mu\text{l}/\text{min}$

**MS Parameters**

Collision Gas (mbar) = 3.24e-3  
Collision Energy (eV) = 5



Reagent

---

**13C3 HFPO-DA\_00009**

**PRODUCT CODE:**

M3HFPO-DA

LOT NUMBER: M3HFPODA0817


COMPOUND:

2,3,3,3-Tetrafluoro-2-(1,1,2,2,3,3,3-heptafluoropropoxy)-<sup>13</sup>C<sub>n</sub>-propanoic acid

**STRUCTURE:**

**CAS #:**

Not available



**MOLECULAR FORMULA:**

 $^{13}\text{C}_3\text{ }^{12}\text{C}_3\text{HF}_{11}\text{O}_3$ 

**CONCENTRATION:**

 $50 \pm 2.5 \mu\text{g/ml}$ 

**CHEMICAL PURITY:**

>98%

**LAST TESTED:** (mm/dd/yyyy)

08/17/2017

**EXPIRY DATE:** (mm/dd/yyyy)

08/17/2020

**RECOMMENDED STORAGE:**

Store ampoule in a cool, dark place

**MOLECULAR WEIGHT:**

333.03

**SOLVENT(S):**

Methanol

**ISOTOPIC PURITY:**

>99%  $^{13}\text{C}$

 $(^{13}\text{C}_3)$ 

~~DOCUMENTATION DATA ATTACHED:~~

Figure 1: LC/MS Data (TIC and Mass Spectrum)


Figure 2: LC/MS/MS Data (Selected MRM Transitions)

**ADDITIONAL INFORMATION:**

- See page 2 for further details.
- Contains ~ 1.5% of two constitutional isomers.
- Product is commercially known as GenX.

**FOR LABORATORY USE ONLY; NOT FOR HUMAN OR DRUG USE**

**Certified By:**

  
B.G. Chittim, General Manager

**Date:** 08/25/2017

(mm/dd/yyyy)

**Wellington Laboratories Inc., 345 Southgate Dr. Guelph ON N1G 3M5 CANADA**  
**519-822-2436 • Fax: 519-822-2849 • [Info@well-labs.com](mailto:Info@well-labs.com)**

#### **INTENDED USE:**

The products prepared by Wellington Laboratories Inc. are for laboratory use only. This certified reference material (CRM) was designed to be used as a standard for the identification and/or quantification of the specific chemical compound it contains.

#### **HAZARDS:**

This product should only be used by qualified personnel familiar with its potential hazards and trained in the handling of hazardous chemicals. Due care should be exercised to prevent unnecessary human contact or ingestion. All procedures should be carried out in a well-functioning fume hood and suitable gloves, eye protection, and clothing should be worn at all times. Waste should be disposed of according to national and regional regulations. Safety Data Sheets (SDSs) are available upon request.

#### **SYNTHESIS / CHARACTERIZATION:**

Where possible, all of our products are synthesized using single-product unambiguous routes. They are then characterized, and their structures and purities confirmed, using a combination of the most relevant techniques, such as NMR, GC/MS, LC/MS/MS, SFC/UV/MS/MS, x-ray crystallography, and melting point. Isotopic purities of mass-labelled compounds are also confirmed using HRGC/HRMS and/or LC/MS/MS.

#### **HOMOGENEITY:**

Prior to solution preparation, crystalline material is tested for homogeneity using a variety of techniques (as stated above) and its solubility in a given diluent is taken into consideration. Duplicate solutions of a new product are prepared from the same crystalline lot and, after the addition of an appropriate internal standard, they are compared by GC/MS, LC/MS/MS and/or SFC/UV/MS/MS. The relative response factors of the analyte of interest in each solution are required to be <5% RSD. New solution lots of existing products are compared to older lots in the same manner, which further confirms the homogeneity of the crystalline material as well as the stability and homogeneity of the solutions in the storage containers. In order to maintain the integrity of the assigned value(s), and associated uncertainty, the dilution or injection of a subsample of this product should be performed using calibrated measuring equipment.

#### **UNCERTAINTY:**

The maximum combined relative standard uncertainty of our reference standard solutions is calculated using the following equation:

The combined relative standard uncertainty,  $u_c(y)$ , of a value  $y$  and the uncertainty of the independent parameters

$x_1, x_2, \dots, x_n$  on which it depends is:

$$u_c(y(x_1, x_2, \dots, x_n)) = \sqrt{\sum_{i=1}^n u(y, x_i)^2}$$

where  $x$  is expressed as a relative standard uncertainty of the individual parameter.

The individual uncertainties taken into account include those associated with weights (calibration of the balance) and volumes (calibration of the volumetric glassware). An expanded maximum combined percent relative uncertainty of  $\pm 5\%$  (calculated with a coverage factor of 2 and a level of confidence of 95%) is stated on the Certificate of Analysis for all of our products.

#### **TRACEABILITY:**

All reference standard solutions are traceable to specific crystalline lots. The microbalances used for solution preparation are regularly tested by an external ISO/IEC 17025 accredited calibration company. In addition, their calibration is verified prior to each weighing using calibrated NIST and/or NRC traceable external weights. All volumetric glassware used is calibrated, of Class A tolerance, and has been tested according to the appropriate ASTM procedures, which are ultimately traceable to NIST. For certain products, traceability to international interlaboratory studies has also been established.

#### **EXPIRY DATE / PERIOD OF VALIDITY:**

Ongoing stability studies of this product have demonstrated stability in its composition and concentration, until the specified expiry date, in the unopened ampoule. Monitoring for any degradation or change in concentration of the listed analyte(s) is performed on a routine basis.

#### **LIMITED WARRANTY:**

At the time of shipment, all products are warranted to be free of defects in material and workmanship and to conform to the stated technical and purity specifications.



#### **QUALITY MANAGEMENT:**

This product was produced using a Quality Management System registered to the latest versions of ISO 9001 by SAI Global, ISO/IEC 17025 by the Canadian Association for Laboratory Accreditation Inc. (CALA; A 1226), and ISO GUIDE 34 by ANSI-ASQ National Accreditation Board (ANAB; AR-1523).



\*\*For additional information or assistance concerning this or any other products from Wellington Laboratories Inc., please visit our website at [www.well-labs.com](http://www.well-labs.com) or contact us directly at [info@well-labs.com](mailto:info@well-labs.com)\*\*

**Figure 1: M3HFPO-DA; LC/MS Data (TIC and Mass Spectrum)**



**Conditions for Figure 1:**

**LC:** Waters Acquity Ultra Performance LC  
**MS:** Micromass Quattro *micro* API MS

**Chromatographic Conditions**

Column: Acquity UPLC BEH Shield RP<sub>18</sub>  
1.7  $\mu$ m, 2.1 x 100 mm

Mobile phase: Gradient  
Start: 55% MeOH / 45% H<sub>2</sub>O with 10 mM NH<sub>4</sub>OAc buffer  
Ramp to 90% organic over 7.5 min and hold for 1.5 min  
before returning to initial conditions in 0.5 min.

Time: 10 min


Flow: 300  $\mu$ l/min

**MS Parameters**

Experiment: Full Scan (225 - 850 amu)

Source: Electrospray (negative)  
Capillary Voltage (kV) = 3.00  
Cone Voltage (V) = 10.00  
Cone Gas Flow (l/hr) = 100  
Desolvation Gas Flow (l/hr) = 750

**Figure 2: M3HFPO-DA; LC/MS/MS Data (Selected MRM Transitions)**



**Conditions for Figure 2:**

Injection: Direct loop Injection  
10  $\mu\text{l}$  (500 ng/ml M3HFPO-DA)

Mobile phase: Isocratic 80% MeOH / 20%  $\text{H}_2\text{O}$  with 10 mM  $\text{NH}_4\text{OAc}$  buffer

Flow: 300  $\mu\text{l}/\text{min}$

**MS Parameters**

Collision Gas (mbar) = 3.24e-3  
Collision Energy (eV) = 5

Reagent

---

**HFPO-DA\_00004**



# WELLINGTON LABORATORIES

## CERTIFICATE OF ANALYSIS DOCUMENTATION

**PRODUCT CODE:**

HFPO-DA

**LOT NUMBER:**


HFPODA0717

**COMPOUND:**

2,3,3,3-Tetrafluoro-2-(1,1,2,2,3,3,3-heptafluoropropoxy)-propanoic acid

**STRUCTURE:****CAS #:**

13252-13-6

**MOLECULAR FORMULA:** $C_6H_2F_{11}O_3$ **MOLECULAR WEIGHT:**

330.05

**CONCENTRATION:** $50 \pm 2.5 \mu\text{g/ml}$ **SOLVENT(S):**

Methanol

**CHEMICAL PURITY:**

&gt;98%

**LAST TESTED:** (mm/dd/yyyy)

07/13/2017

**EXPIRY DATE:** (mm/dd/yyyy)

07/13/2020

**RECOMMENDED STORAGE:**

Store ampoule in a cool, dark place

**DOCUMENTATION/ DATA ATTACHED:**

Figure 1: LC/MS Data (TIC and Mass Spectrum)

Figure 2: LC/MS/MS Data (Selected MRM Transitions)

**ADDITIONAL INFORMATION:**

- See page 2 for further details.
- Product is commercially known as GenX.

**FOR LABORATORY USE ONLY: NOT FOR HUMAN OR DRUG USE**

Certified By:

  
B.G. Chittim, General Manager

Date:

07/14/2017  
(mm/dd/yyyy)

Wellington Laboratories Inc., 345 Southgate Dr. Guelph ON N1G 3M5 CANADA  
519-822-2436 • Fax: 519-822-2849 • info@well-labs.com

#### **INTENDED USE:**

The products prepared by Wellington Laboratories Inc. are for laboratory use only. This certified reference material (CRM) was designed to be used as a standard for the identification and/or quantification of the specific chemical compound it contains.

#### **HAZARDS:**

This product should only be used by qualified personnel familiar with its potential hazards and trained in the handling of hazardous chemicals. Due care should be exercised to prevent unnecessary human contact or ingestion. All procedures should be carried out in a well-functioning fume hood and suitable gloves, eye protection, and clothing should be worn at all times. Waste should be disposed of according to national and regional regulations. Safety Data Sheets (SDSs) are available upon request.

#### **SYNTHESIS / CHARACTERIZATION:**

Where possible, all of our products are synthesized using single-product unambiguous routes. They are then characterized, and their structures and purities confirmed, using a combination of the most relevant techniques, such as NMR, GC/MS, LC/MS/MS, SFC/UV/MS/MS, x-ray crystallography, and melting point. Isotopic purities of mass-labelled compounds are also confirmed using HRGC/HRMS and/or LC/MS/MS.

#### **HOMOGENEITY:**

Prior to solution preparation, crystalline material is tested for homogeneity using a variety of techniques (as stated above) and its solubility in a given diluent is taken into consideration. Duplicate solutions of a new product are prepared from the same crystalline lot and, after the addition of an appropriate internal standard, they are compared by GC/MS, LC/MS/MS and/or SFC/UV/MS/MS. The relative response factors of the analyte of interest in each solution are required to be <5% RSD. New solution lots of existing products are compared to older lots in the same manner, which further confirms the homogeneity of the crystalline material as well as the stability and homogeneity of the solutions in the storage containers. In order to maintain the integrity of the assigned value(s), and associated uncertainty, the dilution or injection of a subsample of this product should be performed using calibrated measuring equipment.

#### **UNCERTAINTY:**

The maximum combined relative standard uncertainty of our reference standard solutions is calculated using the following equation:

The combined relative standard uncertainty,  $u_c(y)$ , of a value  $y$  and the uncertainty of the independent parameters

$x_1, x_2, \dots, x_n$  on which it depends is:

$$u_c(y(x_1, x_2, \dots, x_n)) = \sqrt{\sum_{i=1}^n u(y, x_i)^2}$$

where  $x$  is expressed as a relative standard uncertainty of the individual parameter.

The individual uncertainties taken into account include those associated with weights (calibration of the balance) and volumes (calibration of the volumetric glassware). An expanded maximum combined percent relative uncertainty of  $\pm 5\%$  (calculated with a coverage factor of 2 and a level of confidence of 95%) is stated on the Certificate of Analysis for all of our products.

#### **TRACEABILITY:**

All reference standard solutions are traceable to specific crystalline lots. The microbalances used for solution preparation are regularly tested by an external ISO/IEC 17025 accredited calibration company. In addition, their calibration is verified prior to each weighing using calibrated NIST and/or NRC traceable external weights. All volumetric glassware used is calibrated, of Class A tolerance, and has been tested according to the appropriate ASTM procedures, which are ultimately traceable to NIST. For certain products, traceability to International Interlaboratory studies has also been established.

#### **EXPIRY DATE / PERIOD OF VALIDITY:**

Ongoing stability studies of this product have demonstrated stability in its composition and concentration, until the specified expiry date, in the unopened ampoule. Monitoring for any degradation or change in concentration of the listed analyte(s) is performed on a routine basis.

#### **LIMITED WARRANTY:**

At the time of shipment, all products are warranted to be free of defects in material and workmanship and to conform to the stated technical and purity specifications.

#### **QUALITY MANAGEMENT:**



This product was produced using a Quality Management System registered to the latest versions of ISO 9001 by SAI Global, ISO/IEC 17025 by the Canadian Association for Laboratory Accreditation Inc. (CALA; A 1226), and ISO GUIDE 34 by ANSI-ASQ National Accreditation Board (ANAB; AR-1523).



\*\*For additional information or assistance concerning this or any other products from Wellington Laboratories Inc., please visit our website at [www.well-labs.com](http://www.well-labs.com) or contact us directly at [info@well-labs.com](mailto:info@well-labs.com)\*\*



**Figure 1:** HFPO-DA; LC/MS Data (TIC and Mass Spectrum)



**Conditions for Figure 1:**

**LC:** Waters Acquity Ultra Performance LC  
**MS:** Micromass Quattro *micro* API MS

**Chromatographic Conditions**

Column: Acquity UPLC BEH Shield RP<sub>18</sub>  
1.7  $\mu$ m, 2.1 x 100 mm

Mobile phase: Gradient

Start: 55% MeOH / 45% H<sub>2</sub>O with 10 mM NH<sub>4</sub>OAc buffer  
Ramp to 90% organic over 7.5 min and hold for 1.5 min  
before returning to initial conditions in 0.5 min.

Time: 10 min

Flow: 300  $\mu$ l/min

**MS Parameters**

Experiment: Full Scan (250 - 850 amu)

Source: Electrospray (negative)


Capillary Voltage (kV) = 3.00

Cone Voltage (V) = 10.00

Cone Gas Flow (l/hr) = 100

Desolvation Gas Flow (l/hr) = 700

**Figure 2: HFPO-DA; LC/MS/MS Data (Selected MRM Transitions)**



**Conditions for Figure 2:**

Injection: Direct loop Injection  
10  $\mu$ l (500 ng/ml HFPO-DA)

Mobile phase: Isocratic 80% MeOH / 20% H<sub>2</sub>O with 10 mM NH<sub>4</sub>OAc buffer

Flow: 300  $\mu$ l/min

**MS Parameters**

Collision Gas (mbar) = 3.20e-3  
Collision Energy (eV) = 5

# 8321A\_HFPO\_Du

---

HFPO-DA