

TX-97-241  
HAEL: 97-0216  
EAN: 972790



MUTAGENICITY TEST WITH  
EC 97-0216, Crude MCHM  
IN THE *SALMONELLA-ESCHERICHIA COLI*/MAMMALIAN-MICROSOME REVERSE  
MUTATION ASSAY WITH A CONFIRMATORY ASSAY

FINAL REPORT

AUTHOR

Timothy E. Lawlor, M.A.

SPONSOR'S REPRESENTATIVE

Eugene D. Barber, Ph. D.

PERFORMING LABORATORY

Covance Laboratories Inc. (Covance)  
9200 Leesburg Pike  
Vienna, Virginia 22182

LABORATORY PROJECT ID

Covance Study No.: 18779-0-409R

SUBMITTED TO

Eastman Chemical Company  
Bldg. 320, Kodak Park  
Rochester, NY 14652-6272

STUDY COMPLETION DATE

September 12, 1997

## QUALITY ASSURANCE STATEMENT


STUDY TITLE: *Salmonella - Escherichia coli*/Mammalian-Microsome Reverse Mutation Assay with a Confirmatory Assay

ASSAY NO.: 18779-0-409R

PROTOCOL NO.: 409R, Edition 4, Modified for Eastman Chemical Company

Quality Assurance inspections of the study and review of the final report of the above referenced project were conducted according to the Standard Operating Procedures of the Quality Assurance Unit and according to the general requirements of the appropriate Good Laboratory Practice regulations. Findings from the inspections and final report review were reported to management and to the study director on the following dates:

<u>Inspection/Date</u>	<u>Findings Reported</u>	<u>Auditor</u>
Dilution of Test Article - 08/12/97	08/12/97	S. Buresh
Draft Report Review - 09/02,03/97	09/03/97	K. Cutchins
Final Report Review - 09/12/97	09/12/97	K. Cutchins



Quality Assurance Unit


9/12/97

Date Released

### STUDY COMPLIANCE AND CERTIFICATION

The described study was conducted in compliance with the Good Laboratory Practice regulations as set forth in the Code of Federal Regulations (21 CFR 58, 40 CFR 792, and 40 CFR 160) and Annex 2 of the OECD Guidelines for Testing of Chemicals (C(81)30 Final) as required by Council Directive 87/18/EEC of December 18, 1986. The study described herein also complies with the EEC Annex V Guideline number B.14, "Other Effects-Mutagenicity, Salmonella Typhimurium-Reverse Mutation Assay", and Guideline number B.13, "Other Effects-Mutagenicity, Escherichia Coli-Reverse Mutation Assay." There were no significant deviations from the aforementioned regulations or the signed protocol that would affect the integrity of the study or the interpretation of the test results. The stability of the test article under the conditions of administration was the responsibility of the Sponsor. The raw data have been reviewed by the Study Director, who certifies that the evaluation of the test article as presented herein represents an appropriate conclusion within the context of the study design and evaluation criteria.

Study Director:

  
\_\_\_\_\_  
Timothy E. Lawlor, M.A.  
Bacterial Mutagenesis  
Genetic and Cellular Toxicology

9.12.97  
Study Completion Date

TABLE OF CONTENTS

	Page No.
I. SUMMARY .....	5
II. STUDY INFORMATION .....	7
III. MATERIALS AND METHODS .....	9
IV. RESULTS AND CONCLUSIONS .....	23
V. DATA TABLES .....	26

**SECTION I. SUMMARY**

**INTRODUCTION AND CONCLUSIONS**

## SUMMARY

### A. Introduction

At the request of Eastman Chemical Company, Covance Laboratories Inc. investigated the test article for mutagenic activity in the *Salmonella* - *Escherichia coli*/Mammalian-Microsome Reverse Mutation Assay with a Confirmatory Assay. This assay evaluated the test article and/or its metabolites for their ability to induce reverse mutations at the histidine locus in the genome of specific *Salmonella typhimurium* tester strains and at the tryptophan locus in an *Escherichia coli* tester strain both in the presence and absence of an exogenous metabolic activation system of mammalian microsomal enzymes derived from Aroclor™-induced rat liver (S9).

The doses tested in the mutagenicity assay were selected based on the results of a dose rangefinding study using tester strains TA100 and WP2uvrA(pKM101) and ten doses of test article ranging from 5,000 to 6.67 µg per plate, one plate per dose, both in the presence and absence of S9 mix.

The tester strains used in the mutagenicity study were *Salmonella typhimurium* tester strains TA98, TA100, TA1535, TA1537, and *Escherichia coli* tester strain WP2uvrA(pKM101). The assay was conducted with six doses of test article in both the presence and absence of S9 mix with concurrent vehicle and positive controls using three plates per dose. The doses tested were 5,000, 2,500, 1,000, 500, 250, and 100 µg per plate in both the presence and absence of S9 mix. The results of the initial mutagenicity assay were confirmed in an independent experiment.

### B. Conclusions

The results of the *Salmonella* - *Escherichia coli*/Mammalian-Microsome Reverse Mutation Assay with a Confirmatory Assay indicate that under the conditions of this study, Eastman Chemical Company's test article did not cause a positive increase in the number of revertants per plate of any of the tester strains either in the presence or absence of microsomal enzymes prepared from Aroclor™-induced rat liver (S9).

**SECTION II. STUDY INFORMATION**

**STUDY INFORMATION**

- A. Sponsor: **Eastman Chemical Company**
- B. Test Article: **EC 97-0216, Crude MCHM  
Lot 6-97**
1. Physical Description: **clear colorless liquid**
2. Date Received: **08/04/97**
- C. Type of Assay: *Salmonella - Escherichia coli*/Mammalian-Microsome Reverse  
Mutation Assay with a Confirmatory Assay
1. Protocol Number: Covance Protocol 409R, Edition 4,  
Modified for Eastman Chemical Company
2. Covance Study Number: **18779-0-409R**
- D. Study Dates
1. Study Initiation Date: **08/05/97**
2. Experimental Start Date: **08/07/97**
3. Experimental Termination Date: **08/20/97**
- E. Study Supervisory Personnel
- |                        |                         |
|------------------------|-------------------------|
| Study Director:        | Timothy E. Lawlor, M.A. |
| Associate Scientist:   | Michael S. Mecchi, B.S. |
| Laboratory Supervisor: | Carlos E. Orantes, B.S. |

NOTE: Effective 02 January 1997, the company name was changed from Corning Hazleton Inc. to Covance Laboratories Inc. Study data and documentation may contain either one or both company names or forms of those names (e.g., acronyms).



### **SECTION III. MATERIALS AND METHODS**

## MATERIALS AND METHODS

The experimental materials, methods and procedures are based on those described by Ames *et al* (1975) and Green and Muriel (1976).

### MATERIALS

#### A. Tester Strains

##### 1. *Salmonella typhimurium*

The tester strains used were the *Salmonella typhimurium* histidine auxotrophs TA98, TA100, TA1535, and TA1537 as described by Ames *et al* (1975). The specific genotypes of these strains are shown in Table I.

**TABLE I. TESTER STRAIN GENOTYPES**

Histidine Mutation			Additional Mutations		
<i>his</i> G46	<i>his</i> C3076	<i>his</i> D3052	LPS	Repair	R Factor
TA1535	TA1537		<i>rfa</i>	<i>uvrB</i>	-
TA100		TA98	<i>rfa</i>	<i>uvrB</i>	+R

In addition to a mutation in the histidine operon, the tester strains contain two additional mutations which enhance their sensitivity to some mutagenic compounds. The *rfa* wall mutation results in the loss of one of the enzymes responsible for the synthesis of part of the lipopolysaccharide barrier that forms the surface of the bacterial cell wall. The resulting cell wall deficiency increases permeability to certain classes of chemicals such as those containing large ring systems (i.e. benzo(a)pyrene) that would otherwise be excluded by a normal intact cell wall.

The second mutation, a deletion of the *uvrB* gene, results in a deficient DNA excision repair system which greatly enhances the sensitivity of these strains to some mutagens. Since the *uvrB* deletion extends through the *bio* gene, all of the tester strains containing this deletion require the vitamin biotin for growth.

Strains TA98 and TA100 also contain the R-factor plasmid, pKM101, which further increases the sensitivity of these strains to some mutagens. The mechanism by which this plasmid increases sensitivity to mutagens has been suggested to be by modifying an existing bacterial DNA repair polymerase complex involved with the mismatch-repair process.

Tester strains TA98 and TA1537 are reverted from histidine dependence (auxotrophy) to histidine independence (prototrophy) by frameshift mutagens. TA1535 is reverted by base

substitution mutagens and TA100 is reverted by mutagens which cause both frameshift and base substitution mutations.

2. *Escherichia coli*

The tester strain used was the tryptophan auxotroph WP2*uvrA* as described by Green and Muriel (1976) containing the pKM101 plasmid.

In addition to a mutation in the tryptophan operon, the tester strain contains a *uvrA* DNA repair deficiency which enhances its sensitivity to some mutagenic compounds. This deficiency allows the strain to show enhanced mutability since the *uvrA* repair system would normally act to remove the damaged part of the DNA molecule and accurately repair it afterwards.

Tester strain WP2*uvrA*(pKM101) also contains the R-factor plasmid, pKM101, which further increases the sensitivity of this strain to some mutagens. The mechanism by which this plasmid increases sensitivity to mutagens has been suggested to be by modifying an existing bacterial DNA repair polymerase complex involved with the mismatch-repair process.

Tester strain WP2*uvrA*(pKM101) is reverted from tryptophan dependence (auxotrophy) to tryptophan independence (prototrophy) by base substitution mutagens.

3. Source of Tester Strains

a. *Salmonella typhimurium*

The tester strains in use at Covance were received directly from Dr. Bruce Ames, Department of Biochemistry, University of California, Berkeley.

b. *Escherichia coli*

The tester strain, WP2*uvrA*(pKM101), in use at Covance was received from the National Collection of Industrial Bacteria, Torrey Research Station, Scotland (United Kingdom).

4. Storage of the Tester Strains

a. Frozen Permanent Stocks

Frozen permanent stocks were prepared by growing fresh overnight cultures, adding DMSO (0.09 ml/ml of culture) and freezing small aliquots (0.5-1.5 ml) at  $\leq -70^{\circ}\text{C}$ .

b. Master Plates

Master plates were prepared by streaking each tester strain from a frozen permanent stock onto minimal agar appropriately supplemented with 1) for *Salmonella typhimurium*, an excess of histidine, and biotin, and for tester strains TA98 and TA100, ampicillin (25 µg/ml), to ensure the stable maintenance of the pKM101 plasmid; and 2) for *Escherichia coli*, an excess of tryptophan. Tester strain master plates were stored at  $5 \pm 3^\circ\text{C}$ .

5. Preparation of Overnight Cultures

a. Inoculation

Overnight cultures for use in all testing procedures were inoculated by transferring a colony from the appropriate master plate to a flask containing culture medium. Inoculated flasks were placed in a shaker/incubator which was programmed to begin operation (shaking,  $125 \pm 25$  rpm; incubation,  $37 \pm 2^\circ\text{C}$ ) so that the overnight cultures were in log phase or late log phase when turbidity monitoring began.

b. Harvest

To ensure that cultures were harvested in late log phase, the length of incubation was determined by spectrophotometric monitoring of culture turbidity. Cultures were harvested once a predetermined turbidity was reached as determined by a percent transmittance (%T) reading on a spectrophotometer. This target turbidity ensures that cultures have reached a density of at least  $0.5 \times 10^9$  cells per ml and that the cultures have not overgrown. Overgrown (stationary) cultures may exhibit decreased sensitivity to some mutagens. Cultures were removed from incubation when the target %T was reached and were placed at  $5 \pm 3^\circ\text{C}$ .

6. Confirmation of Tester Strain Genotypes

Tester strain cultures were checked for the following genetic markers on the day of their use in the mutagenicity assay:

a. *Salmonella typhimurium*

1) *rfa* Wall Mutation

The presence of the *rfa* wall mutation was confirmed by demonstration of the sensitivity of the culture to crystal violet. An aliquot of an overnight culture of each strain was overlaid onto plates containing selective media and an antibiotic

sensitivity disk containing 10 µg of crystal violet was added. Sensitivity was demonstrated by inhibition of bacterial growth in a zone immediately surrounding the disk.

2) pKM101 Plasmid R-factor

The presence of the pKM101 plasmid was confirmed for tester strains TA98 and TA100 by demonstration of resistance to ampicillin. An aliquot of an overnight culture of each strain was overlaid onto plates containing selective media and an antibiotic sensitivity disk containing 10 µg of ampicillin was added. Resistance was demonstrated by bacterial growth in the zone immediately surrounding the disk.

3) Characteristic Number of Spontaneous Revertants

The mean number of spontaneous revertants per plate in the vehicle controls that are characteristic of the respective strains were demonstrated by plating 100 µl aliquots of the culture along with the appropriate vehicle on selective media.

b. *Escherichia coli*

1) pKM101 Plasmid R-factor

The presence of the pKM101 plasmid was confirmed for tester strain WP2uvrA(pKM101) by demonstration of resistance to ampicillin. An aliquot of an overnight culture was overlaid onto a plate containing selective media and an antibiotic sensitivity disk containing 10 µg of ampicillin was added. Resistance was demonstrated by bacterial growth in the zone surrounding the disk.

2) Characteristic Number of Spontaneous Revertants

The mean number of spontaneous revertants per plate in the vehicle controls that is characteristic of the strain was demonstrated by plating 100 µl aliquots of the WP2uvrA(pKM101) culture along with the appropriate vehicle on selective media.

7. Tester Strain Media

a. Culturing Broth

The broth used to grow overnight cultures of the tester strains was Vogel-Bonner salt solution (Vogel and Bonner, 1956) supplemented with 2.5% (w/v) Oxoid Nutrient Broth No. 2 (dry powder).

## b. Agar Plates

Bottom agar (25 ml per 15 x 100 mm petri dish) was Vogel-Bonner minimal medium E (Vogel and Bonner, 1956), supplemented with 1.5% (w/v) agar and 0.2% (w/v) glucose.

## c. Overlay Agar for Selection of Revertants

Overlay (top) agar was prepared with 0.7% agar (w/v) and 0.5% NaCl (w/v) and was supplemented with 10 ml of 1) 0.5 mM histidine/biotin solution per 100 ml agar for selection of histidine revertants, or 2) 0.5 mM tryptophan solution per 100 ml of agar for selection of tryptophan revertants. When S9 mix was required, 2.0 ml of the supplemented top agar was used in the overlay. However, when S9 mix was not required, water was added to the supplemented top agar (0.5 ml of water per 2 ml of supplemented top agar) and the resulting 2.5 ml of diluted supplemented top agar was used for the overlay. This dilution ensured that the final top agar and amino acid supplement concentrations remained the same both in the presence and absence of S9 mix.

**B. Liver Microsomal Enzyme Reaction Mixture (S9 Mix)**1. S9 Homogenate

Liver microsomal enzymes (S9 homogenate) were purchased from Molecular Toxicology, Inc., Batch 0755 (35.03 mg of protein per ml). The homogenate was prepared from male Sprague-Dawley rats that had been injected (i.p.) with Aroclor™ 1254 (200 mg per ml in corn oil) at 500 mg/kg as described by Ames *et al*, 1975.

2. S9 Mix

The S9 mix was prepared immediately prior to its use in any experimental procedure. The S9 mix contained the components indicated in Table II.

**TABLE II. S9 MIX COMPONENTS**

H <sub>2</sub> O	0.70 ml
1M NaH <sub>2</sub> PO <sub>4</sub> /Na <sub>2</sub> HPO <sub>4</sub> , pH 7.4	0.10 ml
0.25M Glucose-6-phosphate	0.02 ml
0.10M NADP	0.04 ml
0.825M KCl/0.2M MgCl <sub>2</sub>	0.04 ml
S9 Homogenate	<u>0.10 ml</u>
	1.00 ml

**C. Controls****1. Vehicle Controls**

Vehicle controls were plated for all tester strains both in the presence and absence of S9 mix. The vehicle control was plated, using a 50 µl aliquot of vehicle (equal to the maximum aliquot of test article dilution plated), along with a 100 µl aliquot of the appropriate tester strain and a 500 µl aliquot of S9 mix (when necessary), on selective agar.

**2. Positive Controls**

The combinations of positive controls, activation condition and tester strains plated concurrently with the assay are indicated in Table III.

**TABLE III. POSITIVE CONTROLS**

<u>Tester Strain</u>	<u>S9 Mix</u>	<u>Positive Control</u>	<u>Conc per plate</u>
TA98	+	2-aminoanthracene	2.5 µg
TA98	-	2-nitrofluorene	1.0 µg
TA100	+	2-aminoanthracene	2.5 µg
TA100	-	sodium azide	2.0 µg
TA1535	+	2-aminoanthracene	2.5 µg
TA1535	-	sodium azide	2.0 µg
TA1537	+	2-aminoanthracene	2.5 µg
TA1537	-	ICR-191	2.0 µg
WP2uvrA(pKM101)	+	2-aminoanthracene	5.0 µg
WP2uvrA(pKM101)	-	4-nitroquinoline-N-oxide	2.0 µg

**a. Source and Grade of Positive Control Articles**

2-aminoanthracene (CAS #613-13-8), Sigma Chemical Co., purity ≥ 97.5%; 2-nitrofluorene (CAS #607-57-8), Aldrich Chemical Co., purity 98%; sodium azide (CAS #26628-22-8), Sigma Chemical Co., purity >98%; ICR-191 (CAS #1707-45-0), Sigma Chemical Co., purity 98%; 4-nitroquinoline-N-oxide (CAS #56-57-5), Sigma Chemical Co., purity >99%.

3. Sterility Controls

a. Test Article

The most concentrated test article dilution was checked for sterility by plating a 50 µl aliquot (the same volume used in the assay) on selective agar.

b. S9 Mix

The S9 mix was checked for sterility by plating 0.5 ml on selective agar.

## METHODS

A. Dose Ranging Study

The growth inhibitory effect (cytotoxicity) of the test article to the test system was determined in order to allow the selection of appropriate doses to be tested in the mutagenicity assay.

1. Design

The dose ranging study was performed using tester strains TA100 and WP2uvrA(pKM101) both in the presence and absence of S9 mix. Ten doses of test article were tested at one plate per dose. The test article was checked for cytotoxicity up to a maximum concentration of 5 mg per plate if solubility/miscibility permitted.

a. Rationale

The cytotoxicity of the test article observed on tester strain TA100 is generally representative of that observed on the other tester strains and because of TA100's comparatively high number of spontaneous revertants per plate, gradations of cytotoxicity can be readily discerned from routine experimental variation. The *Escherichia coli* tester strain WP2uvrA(pKM101) does not possess the *rfa* wall mutation that the *Salmonella typhimurium* strains have and thus, a different range of cytotoxicity may be observed. Also, the cytotoxicity induced by a test article in the presence of S9 mix may vary greatly from that observed in the absence of S9 mix. Therefore, this would require that different test article dose ranges be tested in the mutagenicity assay based on the presence or absence of the S9 mix.



2. Evaluation of the Dose Ranging Study

Cytotoxicity is detectable as a decrease in the number of revertant colonies per plate and/or by a thinning or disappearance of the bacterial background lawn.

3. Selection of the Maximum Dose for the Mutagenicity Assay

a. Cytotoxicity Observed

Cytotoxicity was observed in the dose ranging study and the highest concentration of test article used in the subsequent mutagenicity assay was a dose which gave a reduction of revertants per plate and/or a thinning of the bacterial background lawn.

**B. Mutagenicity Assay**

1. Design

The assay was performed using tester strains TA98, TA100, TA1535, TA1537, and WP2uvrA(pKM101) both in the presence and absence of S9 mix. Six doses of test article were tested along with the appropriate vehicle and positive controls. The doses of test article were selected based on the results of the dose ranging study. The results of the initial mutagenicity assay were confirmed in an independent experiment.

2. Frequency and Route of Administration

The tester strains were exposed to the test article via the plate incorporation methodology originally described by Ames *et al* (1975) and Maron and Ames (1983). This methodology has been shown to detect a wide range of classes of chemical mutagens. In the plate incorporation methodology, the test article, the tester strain and the S9 mix (where appropriate) were combined in molten agar which was overlaid onto a minimal agar plate. Following incubation at  $37 \pm 2^\circ\text{C}$  for  $48 \pm 8$  hr, revertant colonies were counted. All doses of the test article, the vehicle controls and the positive controls were plated in triplicate.

**C. Plating Procedures**

These procedures were used in both the dose ranging study and the mutagenicity assay.

Each plate was labeled with a code which identified the test article, test phase, tester strain, activation condition and dose. The S9 mix and dilutions of the test article were prepared immediately prior to their use.

When S9 mix was not required, 100 µl of tester strain and 50 µl of vehicle or test article dose was added to 2.5 ml of molten selective top agar (maintained at  $45 \pm 2^\circ\text{C}$ ). When S9 mix was required, 500 µl of S9 mix, 100 µl of tester strain and 50 µl of vehicle or test article dose was added to 2.0 ml of molten selective top agar. After the required components had been added, the mixture was vortexed and overlaid onto the surface of 25 ml of minimal bottom agar contained in a 15 x 100 mm petri dish. After the overlay had solidified, the plates were inverted and incubated for  $48 \pm 8$  hr at  $37 \pm 2^\circ\text{C}$ . Positive control articles were plated using a 50 µl plating aliquot.

#### **D. Scoring the Plates**

Plates which were not evaluated immediately following the incubation period were held at  $5 \pm 3^\circ\text{C}$  until such time that colony counting and bacterial background lawn evaluation could take place.

##### **1. Bacterial Background Lawn Evaluation**

The condition of the bacterial background lawn was evaluated for evidence of cytotoxicity and test article precipitate. Evidence of cytotoxicity was scored relative to the vehicle control plate and was recorded along with the revertant counts for all plates at that dose on the data tables using the code system presented at the end of the Materials and Methods Section.

##### **2. Counting Revertant Colonies**

The number of revertant colonies per plate for the vehicle controls and all plates containing test article were counted manually. The number of revertant colonies per plate for the positive controls were counted by automated colony counter.

#### **E. Analysis of Data**

For all replicate platings, the mean revertants per plate and the standard deviation were calculated. The results of these calculations are presented in tabular form in the Data Tables Section of this report.

### **EVALUATION OF TEST RESULTS**

Before assay data were evaluated, the criteria for a valid assay had to be met.

**A. Criteria For A Valid Assay**

The following criteria were used to determine a valid assay:

1. Tester Strain Integrity: *Salmonella typhimurium*

a. *rfa* Wall Mutation

To demonstrate the presence of the *rfa* wall mutation, tester strain cultures exhibited sensitivity to crystal violet.

b. pKM101 Plasmid

To demonstrate the presence of the R-factor plasmid, pKM101, cultures of tester strains TA98 and TA100 exhibited resistance to ampicillin.

c. Characteristic Number of Spontaneous Revertants

To demonstrate the requirement for histidine, the tester strain cultures exhibited a characteristic number of spontaneous revertants per plate when plated along with the vehicle under selective conditions. The acceptable ranges for the mean vehicle controls were as follows:

TA98	8 - 60
TA100	60 - 240
TA1535	4 - 45
TA1537	2 - 25

2. Tester Strain Integrity : *Escherichia coli*

a. pKM101 Plasmid

To demonstrate the presence of the R-factor plasmid, pKM101, cultures of tester strain WP2uvrA(pKM101) exhibited resistance to ampicillin.

b. Characteristic Number of Spontaneous Revertants

To demonstrate the requirement for tryptophan, the tester strain culture exhibited a characteristic number of spontaneous revertants per plate when plated along with the vehicle under selective conditions. The acceptable range for the WP2uvrA(pKM101) mean vehicle controls was 80 to 350 revertants per plate.

3. Tester Strain Culture Density

To demonstrate that appropriate numbers of bacteria are plated, the density of tester strain cultures were greater than or equal to  $0.5 \times 10^9$  bacteria per ml and/or had reached a target level of turbidity demonstrated to produce cultures with a density greater than or equal to  $0.5 \times 10^9$  bacteria per ml.

4. Positive Control Values

a. Positive Control Values in the Absence of S9 Mix

To demonstrate that the tester strains were capable of identifying a mutagen, the mean value of a positive control for a respective tester strain exhibited at least a 3-fold increase over the mean value of the vehicle control for that strain.

b. Positive Control Values in the Presence of S9 Mix  
(S9 Mix Integrity)

To demonstrate that the S9 mix was capable of metabolizing a promutagen to its mutagenic form(s), the mean value of the positive control for a respective tester strain in the presence of the S9 mix exhibited at least a 3-fold increase over the mean value of the vehicle control for that strain.

An acceptable positive control in the presence of S9 mix for a specific strain was evaluated as having demonstrated both the integrity of the S9 mix and the ability of the tester strain to detect a mutagen.

5. Cytotoxicity

A minimum of three non-toxic doses were required to evaluate assay data.

**B. Criteria For A Positive Response**

Once the criteria for a valid assay had been met, responses observed in the assay were evaluated as follows:

1. Tester Strains TA98, TA100, and WP2uvrA(pKM101)

For a test article to be considered positive, it had to produce at least a 2-fold increase in the mean revertants per plate of at least one of these tester strains over the mean revertants per plate of the appropriate vehicle control. This increase in the mean number of

revertants per plate had to be accompanied by a dose response to increasing concentrations of the test article.

2. Tester Strains TA1535 and TA1537

For a test article to be considered positive, it had to produce at least a 3-fold increase in the mean revertants per plate of at least one of these tester strains over the mean revertants per plate of the appropriate vehicle control. This increase in the mean number of revertants per plate had to be accompanied by a dose response to increasing concentrations of the test article.

## RECORDS TO BE MAINTAINED

All raw data, documentation, records, the protocol, and the final report generated as a result of this study will be archived in the storage facilities of Covance-Vienna for at least one year following submission of the final report to the Sponsor. After the one year period, the Sponsor may elect to have the aforementioned materials retained in the storage facilities of Covance-Vienna for an additional period of time or sent to a storage facility designated by the Sponsor.

## REFERENCES

- Ames, B.N., J. McCann, and E. Yamasaki. Methods for detecting carcinogens and mutagens with the *Salmonella*/Mammalian-Microsome Mutagenicity Test. *Mutation Research* 31:347-364 (1975).
- Brusick, D.J., V.F. Simmon, H.S. Rosenkranz, V.A. Ray, and R.S. Stafford.  
An evaluation of the *Escherichia coli* WP2 and WP2uvrA reverse mutation assay.  
*Mutation Research* 76:169-190 (1980).
- Green, M.H.L. and W.J. Muriel. Mutagen testing using *trp*<sup>+</sup> reversion in *Escherichia coli*.  
*Mutation Research* 38:3-32 (1976).
- Maron, D.M., and B. Ames. Revised methods for the *Salmonella* Mutagenicity Test.  
*Mutation Research* 113:173-215 (1983).
- Vogel, H.J., and D.M. Bonner. Acetylornithinase of *E. coli*: Partial purification and some properties. *J. Biol. Chem.* 218:97-106 (1956).

**BACTERIAL BACKGROUND LAWN EVALUATION CODE**

The condition of the background bacterial lawn is evaluated both macroscopically and microscopically (using a dissecting microscope) for indications of cytotoxicity and test article precipitate as follows:

<u>CODE</u>	<u>DEFINITION</u>	<u>CHARACTERISTICS OF BACKGROUND LAWN</u>
1	Normal	A healthy microcolony lawn.
2	Slightly Reduced	A noticeable thinning of the microcolony lawn and an increase in the size of the microcolonies compared to the vehicle control plate.
3	Moderately Reduced	A marked thinning of the microcolony lawn and an increase in the size of the microcolonies compared to the vehicle control plate.
4	Extremely Reduced	An extreme thinning of the microcolony lawn and an increase in the size of the microcolonies compared to the vehicle control plate.
5	Absent	A complete lack of any microcolony lawn.
6	Obscured by Precipitate	The background bacterial lawn cannot be accurately evaluated due to microscopic and/or macroscopic test article precipitate.

Evidence of macroscopic test article precipitate on the plates is recorded by addition of the following precipitate code to the code number used to evaluate the condition of the background bacterial lawn.

sp	Slight Precipitate	Noticeable macroscopic precipitate on the plate, however, the precipitate does not influence automated counting of the plate.
mp	Moderate Precipitate	The amount of macroscopic precipitate on the plate would interfere with automated counting, thus requiring the plate to be hand counted.
hp	Heavy Precipitate	The large amount of macroscopic precipitate on the plate makes the required hand counting difficult.

Example: 4mp would indicate a plate observed to have an extremely reduced background lawn which had to be counted manually due to the marked amount of macroscopic test article precipitate.

#### **SECTION IV. RESULTS AND CONCLUSIONS**

## RESULTS

### A. Test Article Handling

The test article was stored at room temperature. Dimethylsulfoxide (DMSO, CAS# 67-68-5, Sigma Chemical Co., Lot 27H0411) was used as the vehicle. At 100 mg per ml, which was the most concentrated stock dilution prepared, the test article formed a clear, colorless solution. The test article remained a solution at all subsequent dilutions prepared.

### B. Dose Rangefinding Study

Doses to be tested in the mutagenicity assay were selected based on the results of the dose rangefinding study conducted on the test article using tester strains TA100 and WP2uvrA(pKM101) in both the presence and absence of S9 mix (one plate per dose). Ten doses of test article, from 5,000 to 6.67 µg per plate, were tested and the results are presented in Tables 1 and 2. These data were generated in Experiment 18779-A1. Cytotoxicity was observed with tester strain TA100 at 3,330 µg per plate and above in both the presence and absence of S9 mix as evidenced by the reduced number of revertants per plate and/or the thinning of the bacterial background lawn. With tester strain WP2uvrA(pKM101), cytotoxicity was observed at 5,000 µg per plate in the presence of S9 mix and at 3,330 µg per plate and above in the absence of S9 mix as evidenced by the reduced number of revertants per plate. No precipitate was observed on the plates at any of the doses tested.

### C. Mutagenicity Assay

The mutagenicity assay results for the test article are presented in Tables 3 through 6. These data were generated in Experiments 18779-B1 and 18779-C1. The data are presented as mean revertants per plate ± standard deviation for each treatment and control group (Tables 4 and 6) and as individual plate counts (Tables 3 and 5).

The results of the dose rangefinding study were used to select six doses to be tested in the mutagenicity assay. The doses tested were 5,000, 2,500, 1,000, 500, 250, and 100 µg per plate in both the presence and absence of S9 mix.

In the initial mutagenicity assay, Experiment 18779-B1 (Tables 3 and 4), and in the confirmatory assay, Experiment 18779-C1 (Tables 5 and 6), all data were acceptable and no positive increases in the number of revertants per plate were observed with any of the tester strains either in the presence or absence of S9 mix.

All criteria for a valid study were met.



## CONCLUSIONS

The results of the *Salmonella* - *Escherichia coli*/Mammalian-Microsome Reverse Mutation Assay with a Confirmatory Assay indicate that under the conditions of this study, Eastman Chemical Company's test article did not cause a positive increase in the number of revertants per plate of any of the tester strains either in the presence or absence of microsomal enzymes prepared from Aroclor™-induced rat liver (S9).

**SECTION V. DATA TABLES**

TABLE 1  
DOSE RANGEFINDING STUDY

TEST ARTICLE ID: EC 97-0216, Crude MCHM

EXPERIMENT ID: 18779-A1

DATE PLATED: 07-Aug-97

VEHICLE: DMSO

DATE COUNTED: 11-Aug-97

$\mu\text{g/PLATE}$	TA100 REVERTANTS PER PLATE			
	WITH S9		WITHOUT S9	
	REVERTANTS PER PLATE	BACKGROUND LAWN EVALUATION*	REVERTANTS PER PLATE	BACKGROUND LAWN EVALUATION*
0.00 (Vehicle) (50 $\mu\text{l}$ )	94	1	84	1
Test Article				
6.67	78	1	89	1
10.0	84	1	94	1
33.3	94	1	87	1
66.7	77	1	77	1
100	106	1	72	1
333	81	1	83	1
667	67	1	74	1
1000	87	1	74	1
3330	57	2	55	2
5000	51	2	36	2

\* Background Lawn Evaluation Codes:

1 = normal	2 = slightly reduced	3 = moderately reduced
4 = extremely reduced	5 = absent	6 = obscured by precipitate
sp = slight precipitate	mp = moderate precipitate (requires hand count)	hp = heavy precipitate (requires hand count)

TABLE 2  
DOSE RANGEFINDING STUDY

TEST ARTICLE ID: EC 97-0216, Crude MCHM

EXPERIMENT ID: 18779-A1

DATE PLATED: 07-Aug-97

VEHICLE: DMSO

DATE COUNTED: 11-Aug-97

µg/PLATE	WP2uvrA(pKM101) REVERTANTS PER PLATE			
	WITH S9		WITHOUT S9	
	REVERTANTS PER PLATE	BACKGROUND LAWN EVALUATION*	REVERTANTS PER PLATE	BACKGROUND LAWN EVALUATION*
0.00 (Vehicle) (50 µl)	132	1	80	1
Test Article				
6.67	147	1	90	1
10.0	119	1	98	1
33.3	134	1	92	1
66.7	143	1	107	1
100	144	1	101	1
333	139	1	105	1
667	121	1	86	1
1000	101	1	84	1
3330	95	1	68	1
5000	60	1	66	1

## \* Background Lawn Evaluation Codes:

1 = normal	2 = slightly reduced	3 = moderately reduced
4 = extremely reduced	5 = absent	6 = obscured by precipitate
sp = slight precipitate	mp = moderate precipitate (requires hand count)	hp = heavy precipitate (requires hand count)

TABLE 3  
MUTAGENICITY ASSAY RESULTS  
INDIVIDUAL PLATE COUNTS

TEST ARTICLE I.D.: EC 97-0216, Crude MCHM

EXPERIMENT I.D.: 18779-B1

DATE PLATED: 12-Aug-97

VEHICLE: DMSO

DATE COUNTED: 14-Aug-97

PLATING ALIQUOT: 50 µl

		REVERTANTS PER PLATE															BACKGROUND	
		TA98			TA100			TA1535			TA1537			WP2uvrA(pKM101)			LAWN*	
		DOSE/PLATE																
MICROSOMES: RAT LIVER			1	2	3	1	2	3	1	2	3	1	2	3	1	2	3	
VEHICLE CONTROL			18	15	8	81	97	102	9	11	10	4	3	9	112	115	110	1
TEST ARTICLE	100 µg		11	21	10	90	109	95	5	14	13	5	11	6	99	120	103	1
	250 µg		8	11	14	104	95	92	12	8	11	10	8	6	113	102	122	1
	500 µg		9	7	9	98	100	86	10	14	8	3	3	5	118	106	125	1
	1000 µg		5	5	6	93	98	82	15	7	8	5	2	3	132	118	92	1
	2500 µg		13	9	5	67	64	62	2	5	4	10	6	3	98	85	69	1
	5000 µg		11	4	10	29	37	45	7	7	9	4	1	2	60	48	47	1 <sup>Ⓢ</sup>
POSITIVE CONTROL**			883	842	982	1007	990	921	76	78	125	259	316	263	1063	1135	1179	1
MICROSOMES: NONE																		
VEHICLE CONTROL			9	9	8	87	71	66	4	8	8	5	1	2	108	125	139	1
TEST ARTICLE	100 µg		9	7	9	77	95	79	5	4	12	4	3	6	103	89	90	1
	250 µg		14	10	6	94	73	66	11	11	8	3	2	3	93	99	113	1
	500 µg		11	11	7	83	81	79	9	7	4	5	3	6	107	94	108	1
	1000 µg		4	7	5	64	82	66	10	5	5	1	3	2	89	76	90	1
	2500 µg		9	4	8	65	69	64	6	3	9	5	1	3	80	88	73	1 <sup>Ⓢ</sup>
	5000 µg		5	1	1	38	32	44	1	4	3	1	0	0	29	36	43	2 <sup>Ⓢ</sup>
POSITIVE CONTROL***			99	83	93	264	409	312	380	308	369	370	335	378	1302	1235	1303	1

** TA98	2-aminoanthracene	2.5 µg/plate	*** TA98	2-nitrofluorene	1.0 µg/plate
TA100	2-aminoanthracene	2.5 µg/plate	TA100	sodium azide	2.0 µg/plate
TA1535	2-aminoanthracene	2.5 µg/plate	TA1535	sodium azide	2.0 µg/plate
TA1537	2-aminoanthracene	2.5 µg/plate	TA1537	ICR-191	2.0 µg/plate
WP2uvrA(pKM101)	2-aminoanthracene	5.0 µg/plate	WP2uvrA(pKM101)	4-nitroquinoline-N-oxide	2.0 µg/plate

\* Background Lawn Evaluation Codes:

1 = normal	2 = slightly reduced	3 = moderately reduced
4 = extremely reduced	5 = absent	6 = obscured by precipitate
sp = slight precipitate	mp = moderate precipitate	hp = heavy precipitate
	(requires hand count)	(requires hand count)

Ⓢ Background bacterial lawns for tester strains TA98 and TA100 were evaluated as slightly reduced (2) at this dose.

Ⓢ Background bacterial lawn for tester strain WP2uvrA(pKM101) was evaluated as normal (1) at this dose.

TABLE 4  
MUTAGENICITY ASSAY RESULTS  
SUMMARY

TEST ARTICLE I.D.: EC 97-0216, Crude MCHM

EXPERIMENT I.D.: 18779-B1

DATE PLATED: 12-Aug-97

VEHICLE: DMSO

DATE COUNTED: 14-Aug-97

PLATING ALIQUOT: 50 µl

MEAN REVERTANTS PER PLATE WITH STANDARD DEVIATION											BACKGROUND
											LAWN*
	DOSE/PLATE	TA98		TA100		TA1535		TA1537		WP2uvrA(pKM101)	
		MEAN	S.D.	MEAN	S.D.	MEAN	S.D.	MEAN	S.D.	MEAN	S.D.
MICROSOMES: RAT LIVER											
VEHICLE CONTROL		14	5	93	11	10	1	5	3	112	3
TEST ARTICLE	100 µg	14	6	98	10	11	5	7	3	107	11
	250 µg	11	3	97	6	10	2	8	2	112	10
	500 µg	8	1	95	8	11	3	4	1	116	10
	1000 µg	5	1	91	8	10	4	3	2	114	20
	2500 µg	9	4	64	3	4	2	6	4	84	15
	5000 µg	8	4	37	8	8	1	2	2	52	7
POSITIVE CONTROL**		902	72	973	46	93	28	279	32	1126	59
MICROSOMES: NONE											
VEHICLE CONTROL		9	1	75	11	7	2	3	2	124	16
TEST ARTICLE	100 µg	8	1	84	10	7	4	4	2	94	8
	250 µg	9	3	78	15	10	2	3	1	102	10
	500 µg	10	3	81	2	7	3	5	2	103	8
	1000 µg	8	3	71	10	7	3	2	1	85	8
	2500 µg	6	2	66	3	6	3	3	2	80	8
	5000 µg	5	3	38	6	3	2	0	1	36	7
POSITIVE CONTROL***		92	8	328	74	352	39	361	23	1280	39

** TA98	2-aminoanthracene	2.5 µg/plate	*** TA98	2-nitrofluorene	1.0 µg/plate
TA100	2-aminoanthracene	2.5 µg/plate	TA100	sodium azide	2.0 µg/plate
TA1535	2-aminoanthracene	2.5 µg/plate	TA1535	sodium azide	2.0 µg/plate
TA1537	2-aminoanthracene	2.5 µg/plate	TA1537	ICR-191	2.0 µg/plate
WP2uvrA(pKM101)	2-aminoanthracene	5.0 µg/plate	WP2uvrA(pKM101)	4-nitroquinoline-N-oxide	2.0 µg/plate

\* Background Lawn Evaluation Codes:

1 = normal	2 = slightly reduced	3 = moderately reduced
4 = extremely reduced	5 = absent	6 = obscured by precipitate
sp = slight precipitate	mp = moderate precipitate	hp = heavy precipitate
	(requires hand count)	(requires hand count)

⊗ Background bacterial lawns for tester strains TA98 and TA100 were evaluated as slightly reduced (2) at this dose.

⊙ Background bacterial lawn for tester strain WP2uvrA(pKM101) was evaluated as normal (1) at this dose.

TABLE 5  
MUTAGENICITY ASSAY RESULTS  
INDIVIDUAL PLATE COUNTS

TEST ARTICLE I.D.: EC 97-0216, Crude MCHM

EXPERIMENT I.D.: 18779-C1

DATE PLATED: 15-Aug-97

VEHICLE: DMSO

DATE COUNTED: 19-Aug-97

PLATING ALIQUOT: 50 µl

		REVERTANTS PER PLATE												BACKGROUND			
		TA98			TA100			TA1535			TA1537			WP2uvrA(pKM101)			LAWN*
DOSE/PLATE		1	2	3	1	2	3	1	2	3	1	2	3	1	2	3	
MICROSOMES: RAT LIVER		1	2	3	1	2	3	1	2	3	1	2	3	1	2	3	
VEHICLE CONTROL		27	38	25	111	113	143	10	12	11	7	5	11	244	248	246	1
TEST ARTICLE	100 µg	10	23	24	121	135	133	11	17	13	11	3	7	260	245	241	1
	250 µg	26	23	19	120	132	121	12	16	9	5	5	4	260	259	260	1
	500 µg	18	21	15	135	109	110	14	16	17	8	10	8	244	279	256	1
	1000 µg	22	21	26	114	102	121	20	17	12	6	12	8	241	256	236	1
	2500 µg	14	19	20	95	74	94	5	6	8	10	8	7	192	218	205	1
	5000 µg	22	16	19	79	61	77	4	5	4	5	8	6	122	172	163	2e
POSITIVE CONTROL**		750	959	757	949	1026	1164	107	126	127	157	149	135	940	924	936	1
MICROSOMES: NONE																	
VEHICLE CONTROL		16	13	17	89	108	105	11	11	14	6	3	7	220	216	220	1
TEST ARTICLE	100 µg	8	14	12	91	100	92	19	19	18	5	11	5	249	220	223	1
	250 µg	10	11	8	119	116	102	15	18	16	7	2	5	261	213	231	1
	500 µg	10	16	18	88	96	116	16	15	15	4	4	7	243	202	230	1
	1000 µg	6	15	10	87	94	96	11	10	6	10	5	3	200	278	246	1
	2500 µg	11	10	10	68	71	86	6	5	8	2	3	2	154	144	211	2e
	5000 µg	5	0	0	66	44	46	0	1	0	0	10	1	133	152	136	3e
POSITIVE CONTROL***		168	187	155	743	701	726	564	598	579	292	270	323	1318	1324	1341	1

\*\* TA98 2-aminoanthracene 2.5 µg/plate  
TA100 2-aminoanthracene 2.5 µg/plate  
TA1535 2-aminoanthracene 2.5 µg/plate  
TA1537 2-aminoanthracene 2.5 µg/plate  
WP2uvrA(pKM101) 2-aminoanthracene 5.0 µg/plate

\*\*\* TA98 2-nitrofluorene 1.0 µg/plate  
TA100 sodium azide 2.0 µg/plate  
TA1535 sodium azide 2.0 µg/plate  
TA1537 ICR-191 2.0 µg/plate  
WP2uvrA(pKM101) 4-nitroquinoline-N-oxide 2.0 µg/plate

\* Background Lawn Evaluation Codes:

1 = normal  
2 = slightly reduced  
3 = moderately reduced  
4 = extremely reduced  
5 = absent  
6 = obscured by precipitate  
sp = slight precipitate  
mp = moderate precipitate  
(requires hand count)

3 = moderately reduced  
6 = obscured by precipitate  
hp = heavy precipitate  
(requires hand count)

⊕ Bacterial background lawn for tester strain WP2uvrA(pKM101) was evaluated as normal (1) at this dose.

TABLE 6  
MUTAGENICITY ASSAY RESULTS  
SUMMARY

TEST ARTICLE I.D.: EC 97-0216, Crude MCHM

EXPERIMENT I.D.: 18779-C1

DATE PLATED: 15-Aug-97

VEHICLE: DMSO

DATE COUNTED: 19-Aug-97

PLATING ALIQUOT: 50 µl

MEAN REVERTANTS PER PLATE WITH STANDARD DEVIATION											BACKGROUND LAWN*
DOSE/PLATE	TA98		TA100		TA1535		TA1537		WP2uvrA(pKM101)		
	MEAN	S.D.	MEAN	S.D.	MEAN	S.D.	MEAN	S.D.	MEAN	S.D.	
MICROSOMES: RAT LIVER											
VEHICLE CONTROL	30	7	122	18	11	1	8	3	246	2	1
TEST ARTICLE											
100 µg	19	8	130	8	14	3	7	4	249	10	1
250 µg	23	4	124	7	12	4	5	1	260	1	1
500 µg	18	3	118	15	16	2	9	1	260	18	1
1000 µg	23	3	112	10	16	4	9	3	244	10	1
2500 µg	18	3	88	12	6	2	8	2	205	13	1
5000 µg	19	3	72	10	4	1	6	2	152	27	2*
POSITIVE CONTROL**	822	119	1046	109	120	11	147	11	933	8	1
MICROSOMES: NONE											
VEHICLE CONTROL	15	2	101	10	12	2	5	2	219	2	1
TEST ARTICLE											
100 µg	11	3	94	5	19	1	7	3	231	16	1
250 µg	10	2	112	9	16	2	5	3	235	24	1
500 µg	15	4	100	14	15	1	5	2	225	21	1
1000 µg	10	5	92	5	9	3	6	4	241	39	1
2500 µg	10	1	75	10	6	2	2	1	170	36	2*
5000 µg	2	3	52	12	0	1	4	6	140	10	3*
POSITIVE CONTROL***	170	16	723	21	580	17	295	27	1328	12	1

\*\* TA98 2-aminoanthracene 2.5 µg/plate  
TA100 2-aminoanthracene 2.5 µg/plate  
TA1535 2-aminoanthracene 2.5 µg/plate  
TA1537 2-aminoanthracene 2.5 µg/plate  
WP2uvrA(pKM101) 2-aminoanthracene 5.0 µg/plate

\*\*\* TA98 2-nitrofluorene 1.0 µg/plate  
TA100 sodium azide 2.0 µg/plate  
TA1535 sodium azide 2.0 µg/plate  
TA1537 ICR-191 2.0 µg/plate  
WP2uvrA(pKM101) 4-nitroquinoline-N-oxide 2.0 µg/plate

\* Background Lawn Evaluation Codes:

1 = normal	2 = slightly reduced	3 = moderately reduced
4 = extremely reduced	5 = absent	6 = obscured by precipitate
sp = slight precipitate	mp = moderate precipitate (requires hand count)	hp = heavy precipitate (requires hand count)

\* Bacterial background lawn for tester strain WP2uvrA(pKM101) was evaluated as normal (1) at this dose.