

4-Methylcyclohexane Methanol Study Using XAD Tubes

	Average Front Conc. µg/L	Average Back Conc. µg/L	Total µg/L	Percent Recovery
4-Methylcyclohexane Methanol (Cis & Trans)				
Sample #1	178	3.93	182	91
Sample#2	181	4.87	186	93
Sample #3	185	2.35	187	94
Sample#4	47.3	3.79	51.1	102
Sample #5	51.5	5.16	56.7	113
Sample#6	48.4	3.98	52.4	105

Note: Samples 1-3 Spiked with 40µL of 10 ppm 4-MCHM Std. and extracted with 2 ml of DCM
 Samples 4-6 Spiked with 10µL of 10 ppm 4-MCHM Std. and extracted with 2 ml of DCM.