CETIFICATION

SDG No:

JC22613

Laboratory:

Accutest, New Jersey

Site:

BMS, Building 5 Area, PR

Matrix:

Groundwater

Humacao, PR

SUMMARY:

Groundwater samples (Table 1) were collected on the BMSMC facility – Building 5 Area. The BMSMC facility is located in Humacao, PR. Samples were taken June 16-20, 2016 and were analyzed in Accutest Laboratory of Dayton, New Jersey for the ABN TCL Special List (1,4-Dioxane and Naphthalene were analyzed following the SIM technique); TCL pesticides list; and for low molecular weight alcohols (LMWA) the results were reported under SDG No.: JC22613. Results were validated using the latest validation guidelines (July, 2015) of the EPA Hazardous Waste Support Section. The analyses performed are shown in Table 1. Individual data review worksheets are enclosed for each target analyte group. The data sample organic data samples summary form shows for analytes results that were qualified.

In summary the results are valid and can be used for decision taking purposes.

Table 1. Samples analyzed and analysis performed

SAMPLE ID	SAMPLE DESCRIPTION	MATRIX	ANALYSIS PERFORMED
JC22613-1	MW-3	Groundwater	ABN TCL special list; 1,-4-dioxane and Naphthalene (SIM); LMWA
JC22613-2	MW-5	Groundwater	ABN TCL special list; 1,-4-dioxane and Naphthalene (SIM); LMWA
JC22613-3	MW-13	Groundwater	ABN TCL special list; 1,-4-dioxane and Naphthalene (SIM); LMWA
JC22613-4	BR-4	Groundwater	ABN TCL special list; 1,-4-dioxane and Naphthalene (SIM); LMWA
JC22613-5	MW-16	Groundwater	ABN TCL special list; 1,-4-dioxane and Naphthalene (SIM); LMWA
JC22613-6	MW-9	Groundwater	ABN TCL special list; 1,-4-dioxane and Naphthalene (SIM); LMWA
JC22613-7	MW-23S	Groundwater	ABN TCL special list; 1,-4-dioxane and Naphthalene (SIM); Pesticides TCL list; LMWA
JC22613-7	MW-21S	Groundwater	ABN TCL special list; 1,-4-dioxane and Naphthalene (SIM); Pesticides TCL list; LMWA

Méndez

Reviewer Name:

Rafael Infante

Chemist License 1888

Signature:

Date:

July 21/2016

Report of Analysis

Page 1 of 3

Client Sample ID:	MW-3
Lab Sample ID:	JC22613-1

Matrix: Method:

Project:

AQ - Ground Water

SW846 8270D SW846 3510C BMSMC, Building 5 Area, PR

Date Sampled: 06/16/16 Date Received: 06/21/16

Q

J

Percent Solids: n/a

	File ID	DF	Analyzed	Ву	Prep Date	Prep Batch	Analytical Batch
Run #1	M125545.D	1	06/24/16	SB	06/23/16	OP95018	EM5319

Run #2

Initial Volume Final Volume Run #1 900 ml 1.0 ml

Run #2

ABN TCL Special List

CAS No.	Compound	Result	RL	MDL	Units
95-57-8	2-Chlorophenol	ND	5.6	0.91	ug/I
59-50-7	4-Chloro-3-methyl phenol	ND	5.6	0.99	ug/l
120-83-2	2,4-Dichlorophenol	ND	2.2	1.4	ug/l
105-67-9	2,4-Dimethylphenol	ND	5.6	2.7	ug/l
51-28-5	2,4-Dinitrophenol	ND	11	1.7	ug/l
534-52-1	4,6-Dinitro-o-cresol	ND	5.6	1.4	ug/l
95-48-7	2-Methylphenol	ND	2.2	0.99	ug/l
	3&4-Methylphenol	ND	2.2	0.98	ug/l
88-75-5	2-Nitrophenol	ND	5.6	1.1	ug/l
100-02-7	4-Nitrophenol	ND	11	1.3	ug/l
87-86-5	Pentachlorophenol	ND	5.6	1.5	ug/l
108-95-2	Phenot	ND	2.2	0.44	ug/l
58-90-2	2,3,4,6-Tetrachlorophenol	ND	5.6	1.6	ug/l
95-95-4	2,4,5-Trichlorophenol	ND	5.6	1.5	ug/l
88-06-2	2,4,6-Trichlorophenol	ND	5.6	1.0	ug/I
83-32-9	Acenaphthene	0.50	1.1	0.21	ug/l
208-96-8	Acenaphthylene	ND	1.1	0.15	ug/l
98-86-2	Acetophenone	ND	2.2	0.23	ug/l
120-12-7	Anthracene	ND	1.1	0.23	ug/l
1912-24-9	Atrazine	ND	2.2	0.50	ug/l
100-52-7	Benzaldehyde	ND	5.6	0.32	ug/l
56-55-3	Benzo(a)anthracene	ND	1.1	0.23	ug/l
50-32-8	Benzo(a) pyrene	ND	1.1	0.24	ug/l
205-99-2	Benzo(b)fluoranthene	ND	1.1	0.23	ug/l
191-24-2	Benzo(g,h,i)perylene	ND	1.1	0.38	ug/l
207-08-9	Benzo(k)fluoranthene	ND	1.1	0.23	ug/l
101-55-3	4-Bromophenyl phenyl ether	ND	2.2	0.45	ug/l
85-68-7	Butyl benzyl phthalate	ND	2.2	0.51	ug/l
92-52-4	1,1'-Biphenyl	ND	1.1	0.24	ug/l
91-58-7	2-Chloronaphthalene	ND	2.2	0.26	ug/l
106-47-8	4-Chloroaniline	ND	5.6	0.38	ug/l
86-74-8	Carbazole	ND	1.1	0.25	ug/l

fael Infante Méndez IC # 1888

ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

Report of Analysis

Client Sample ID: Lab Sample ID:	MW-3
Lab Sample ID:	JC22613-1

Matrix: Method: AQ - Ground Water

SW846 8270D SW846 3510C Project: BMSMC, Building 5 Area, PR

Date Sampled: 06/16/16 Date Received: 06/21/16

Percent Solids:

ABN TCL Special List

	•					
CAS No.	Compound	Result	RL	MDL	Units	Q
105-60-2	Caprolactam	ND	2.2	0.72	ug/l	
218-01-9	Chrysene	ND	1.1	0.20	ug/l	
111-91-1	bis(2-Chloroethoxy)methane	ND	2.2	0.31	ug/l	
111-44-4	bis(2-Chloroethyl)ether	ND	2.2	0.28	ug/l	
108-60-1	bis(2-Chloroisopropyl)ether	ND	2.2	0.45	ug/l	
7005-72-3	4-Chlorophenyl phenyl ether	ND	2.2	0.41	ug/l	
121-14-2	2,4-Dinitrotoluene	ND	1.1	0.61	ug/l	
606-20-2	2,6-Dinitrotoluene	ND	1.1	0.53	ug/l	
91-94-1	3,3'-Dichlorobenzidine	ND	2.2	0.56	ug/l	
53-70-3	Dibenzo(a,h)anthracene	ND	1.1	0.37	ug/l	
132-64-9	Dibenzofuran	ND	5.6	0.24	ug/l	
84-74-2	Di-n-butyl phthalate	ND	2.2	0.55	ug/l	
117-84-0	Di-n-octyl phthalate	ND	2.2	0.26	ug/l	
84-66-2	Diethyl phthalate	ND	2.2	0.29	ug/l	
131-11-3	Dimethyl phthalate	ND	2.2	0.24	ug/l	
117-81-7	bis(2-Ethylhexyl)phthalate	ND	2.2	1.8	ug/l	
206-44-0	Fluoranthene	ND	1.1	0.19	ug/l	
86-73-7	Fluorene	1.1	1.1	0.19	ug/l	
118-74-1	Hexachlorobenzene	ND	1.1	0.36	ug/l	
87-68-3	Hexachlorobutadiene	ND	1.1	0.55	ug/l	
77-47-4	Hexachlorocyclopentadiene	ND	11	3.1	ug/l	
67-72-1	Hexachloroethane	ND	2.2	0.43	ug/l	
193-39-5	Indeno(1,2,3-cd)pyrene	ND	1.1	0.37	ug/l	
78-59-1	Isophorone	ND	2.2	0.31	ug/l	
90-12-0	1-Methylnaphthalene	35.9	1.1	0.29	ug/l	
91-57-6	2-Methylnaphthalene	33.3	1.1	0.23	ug/l	
88-74-4	2-Nitroaniline	ND	5.6	0.31	ug/l	
99-09-2	3-Nitroaniline	ND	5.6	0.43	ug/l	
100-01-6	4-Nitroaniline	ND	5.6	0.49	ug/l	
98-95-3	Nitrobenzene	ND	2.2	0.71	ug/l	
621-64-7	N-Nitroso-di-n-propylamine	ND	2.2	0.53	ug/l	
86-30-6	N-Nitrosodiphenylamine	ND	5.6	0.25	ug/l	
85-01-8	Phenanthrene	ND	1.1	0.19	ug/l	
129-00-0	Pyrene	ND	1.1	0.24	ug/l	- 1-
95-94-3	1,2,4,5-Tetrachlorobenzene	ND	2.2	0.41	ug/l	1
CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Lim	its	\
367-12-4	2-Fluorophenol	45%		14-8	8%	
4165-62-2	Phenol-d5	30%		10-1	10%	



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

N = Indicates presumptive evidence of a compound

Report of Analysis

Client Sample ID: MW-3 Lab Sample ID:

JC22613-1

AQ - Ground Water SW846 8270D SW846 3510C Date Sampled: Date Received: 06/21/16 Percent Solids: n/a

06/16/16

Method: Project:

Matrix:

BMSMC, Building 5 Area, PR

ABN TCL Special List

CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limits
118-79-6 4165-60-0 321-60-8	2,4,6-Tribromophenol Nitrobenzene-d5 2-Fluorobiphenyl	97% 72% 74%		39-149% 32-128% 35-119%
1718-51-0	Terphenyl-d14	85%		10-126%



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

321-60-8

1718-51-0

Report of Analysis

Page 1 of 1

Client Sam Lab Sampl Matrix: Method: Project:	le ID: JC226 AQ - C SW846	13-1 Ground Wat 5 8270D BY		3510C		Date		5/16/16 5/21/16 ⁄a
Run #1 Run #2	File ID 3M62544.D	DF 1	Analyzed 06/26/16	By AC	Prep D 06/23/1		Prep Batch OP95018A	Analytical Batch E3M2951
Run #1 Run #2	Initial Volume 900 ml	Final Vo	dune					
CAS No.	Compound		Result	RL	MDL	Units	Q	
91-20-3 123-91-1	Naphthalene 1,4-Dioxane		2.19 ND	0.11 0.11	0.033 0.054	ug/l ug/l		
CAS No.	Surrogate Re	coveries	Run# 1	Run# 2	Lim	its		
4165-60-0	Nitrobenzene-	d5	67%		24-1	25%		

57%

81%



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

2-Fluorobiphenyl

Terphenyl-d14

J = Indicates an estimated value

19-127%

10-119%

B = Indicates analyte found in associated method blank



Report of Analysis

Page 1 of 1

Client Sample ID: MW-3 Lab Sample ID: JC22613-1

Matrix: Method:

Project:

AQ - Ground Water SW846-8015C (DAI)

BMSMC, Building 5 Area, PR

Date Sampled: Date Received:

06/16/16 06/21/16

Percent Solids: n/a

	nalytical Batch GH5333
--	---------------------------

Low Molecular Alcohol List

CAS No.	Compound	Result	RL	MDL	Units	Q
64-17-5	Ethanol	ND	100	55	ug/l	
78-83-1	Isobutyl Alcohol	ND	100	36	ug/l	
67-63-0	Isopropyl Alcohol	ND	100	68	ug/l	
71-23-8	n-Propyl Alcohol	ND	100	43	ug/l	
71-36-3	n-Butyl Alcohol	ND	100	87	ug/l	
78-92-2	sec-Butyl Alcohol	ND	100	66	ug/l	
67-56-1	Methanol	ND	200	71	ug/l	
CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limi	its	
111-27-3	Hexanol	78%		56-1-	45%	





MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

N = Indicates presumptive evidence of a compound

Report of Analysis

SB

06/24/16

Page 1 of 3

Client Sample ID:	MW-5
Client Sample ID: Lab Sample ID:	JC22613-2

Matrix:

AQ - Ground Water

1

Date Sampled: 06/17/16 Date Received: 06/21/16

Method: Project:

SW846 8270D SW846 3510C BMSMC, Building 5 Area, PR

Percent Solids:

Q

File ID Ву DF Analyzed Prep Date

07/01/16

Analytical Batch Prep Batch OP95059 E3E3716

Run #1 Run #2

> Final Volume **Initial Volume**

3E84767.D

Run #1 Run #2 900 ml

1.0 ml

ABN TCL Special List

CAS No.	Compound	Result	RL	MDL	Unit
95-57-8	2-Chlorophenol	ND	5.6	0.91	ug/l
59-50-7	4-Chloro-3-methyl phenol	ND	5.6	0.99	ug/I
120-83-2	2,4-Dichlorophenol	ND	2.2	1.4	ug/l
105-67-9	2,4-Dimethylphenol	ND	5.6	2.7	ug/I
51-28-5	2,4-Dinitrophenol	ND	11	1.7	ug/l
534-52-1	4,6-Dinitro-o-cresol	ND	5.6	1.4	ug/l
95-48-7	2-Methylphenol	ND	2.2	0.99	ug/l
	3&4-Methylphenol	ND	2.2	0.98	ug/l
88-75-5	2-Nitrophenol	ND	5.6	1.1	ug/l
100-02-7	4-Nitrophenol	ND	11	1.3	ug/l
87-86-5	Pentachlorophenol	ND	5.6	1.5	ug/l
108-95-2	Phenol	ND	2.2	0.44	ug/l
58-90-2	2,3,4,6-Tetrachlorophenol	ND	5.6	1.6	ug/i
95-95-4	2,4,5-Trichlorophenol	ND	5.6	1.5	ug/l
88-06-2	2,4,6-Trichlorophenol	ND	5.6	1.0	ug/i
83-32-9	Acenaphthene	ND	1.1	0.21	ug/l
208-96-8	Acenaphthylene	ND	1.1	0.15	ug/l
98-86-2	Acetophenone	ND	2.2	0.23	ug/l
120-12-7	Anthracene	ND	1.1	0.23	ug/l
1912-24-9	Atrazine	ND	2.2	0.50	ug/I
100-52-7	Benzaldehyde	ND	5.6	0.32	ug/l
56-55-3	Benzo(a)anthracene	ND	1.1	0.23	ug/l
50-32-8	Benzo(a) pyrene	ND	1.1	0.24	ug/l
205-99-2	Benzo(b) fluoranthene	ND	1.1	0.23	ug/I
191-24-2	Benzo(g,h,i)perylene	ND	1.1	0.38	ug/l
207-08-9	Benzo(k)fluoranthene	ND	1.1	0.23	ug/l
101-55-3	4-Bromophenyl phenyl ether	ND	2.2	0.45	ug/l
85-68-7	Butyl benzyl phthalate	ND	2.2	0.51	ug/l
92-52-4	1,1'-Biphenyl	ND	1.1	0.24	ug/l
91-58-7	2-Chloronaphthalene	ND	2.2	0.26	ug/l
106-47-8	4-Chloroaniline	ND	5.6	0.38	ug/l
86-74-8	Carbazole	ND	1.1	0.25	ug/l
					-



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

Method:

Project:

Report of Analysis

Client Sample ID: MW-5
Lab Sample ID: JC22613-2
Matrix: AQ - Grou

AQ - Ground Water

SW846 8270D SW846 3510C BMSMC, Building 5 Area, PR

Date Sampled: 06/17/16 Date Received: 06/21/16

Percent Solids: n/a

ABN TCL Special List

ABN ICL	Special List					
CAS No.	Compound	Result	RL	MDL	Units	Q
105-60-2	Caprolactam	ND	2.2	0.72	ug/l	
218-01-9	Chrysene	ND	1.1	0.20	ug/l	
111-91-1	bis(2-Chloroethoxy)methane	ND	2.2	0.31	ug/l	
111-44-4	bis(2-Chloroethyl)ether	ND	2.2	0.28	ug/l	
108-60-1	bis (2-Chloroisopropyl) ether	ND	2.2	0.45	ug/l	
7005-72-3	4-Chlorophenyl phenyl ether	ND	2.2	0.41	ug/l	
121-14-2	2,4-Dinitrotoluene	ND	1.1	0.61	ug/l	
606-20-2	2,6-Dinitrotoluene	ND	1.1	0.53	ug/l	
91-94-1	3,3'-Dichlorobenzidine	ND	2.2	0.56	ug/l	
53-70-3	Dibenzo(a,h)anthracene	ND	1.1	0.37	ug/l	
132-64-9	Dibenzofuran	ND	5.6	0.24	ug/l	
84-74-2	Di-n-butyl phthalate	ND	2.2	0.55	ug/l	
117-84-0	Di-n-octyl phthalate	ND	2.2	0.26	ug/l	
84-66-2	Diethyl phthalate	ND	2.2	0.29	ug/l	
131-11-3	Dimethyl phthalate	ND	2.2	0.24	ug/l	
117-81-7	bis (2-Ethylhexyl) phthalate	ND	2.2	1.8	ug/l	
206-44-0	Fluoranthene	ND	1.1	0.19	ug/l	
86-73-7	Fluorene	ND	1.1	0.19	ug/l	
118-74-1	Hexachlorobenzene	ND	1,1	0.36	ug/l	
87-68-3	Hexachlorobutadiene	ND	1.1	0.55	ug/l	
77-47-4	Hexachlorocyclopentadiene	ND	11	3.1	ug/l	
67-72-1	Hexachloroethane	ND	2.2	0.43	ug/l	
193-39-5	Indeno(1,2,3-cd)pyrene	ND	1.1	0.37	ug/l	
78-59-1	Isophorone	ND	2.2	0.31	ug/l	
90-12-0	1-Methylnaphthalene	ND	1.1	0.29	ug/l	
91-57-6	2-Methylnaphthalene	ND	1.1	0.23	ug/l	
88-74-4	2-Nitroaniline	ND	5.6	0.31	ug/l	
99-09-2	3-Nitroaniline	ND	5.6	0.43	ug/l	
100-01-6	4-Nitroaniline	ND	5.6	0.49	ug/l	
98-95-3	Nitrobenzene	ND	2.2	0.71	ug/l	200
621-64-7	N-Nitroso-di-n-propylamine	ND	2.2	0.53	ug/l	SE NOCKOO DE A
86-30-6	N-Nitrosodiphenylamine	ND	5.6	0.25	ug/l	age
85-01-8	Phenanthrene	ND	1.1	0.19	ug/l	
129-00-0	Pyrene	ND	1.1	0.24	ug/l	Tael Infante
95-94-3	1,2,4,5-Tetrachlorobenzene	ND	2.2	0.41	ug/I	Méndez 10 # 1888
CAS No.	Surrogate Recoveries	Run#1	Run# 2	Lim	its	CO LICENCIA
367-12-4	2-Fluorophenol	62%		14-8	8%	ALCO FICENCIA
4165-62-2	Phenol-d5	40%		10-1	10%	

ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

N = Indicates presumptive evidence of a compound

ACCUTEST

Report of Analysis

Page 3 of 3

Client Sample ID: MW-5

Lab Sample ID: JC22613-2

Matrix:

AQ - Ground Water

Date Sampled: Date Received:

06/17/16 06/21/16

Method: Project: SW846 8270D SW846 3510C BMSMC, Building 5 Area, PR

Percent Solids: n/a

ABN TCL Special List

CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limits
118-79-6	2,4,6-Tribromophenol	67%		39-149%
4165-60-0	Nitrobenzene-d5	94%		32-128%
321-60-8	2-Fluorobiphenyl	89%		35-119%
1718-51-0	Terphenyl-d14	100%		10-126%



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

Report of Analysis

Page 1 of 1

Client Sample ID:	MW-5
Lab Sample ID:	JC22613-2

Matrix:

AQ - Ground Water

Date Sampled: 06/17/16 Date Received: 06/21/16

Method: SW846 8270D BY SIM SW846 3510C Project:

BMSMC, Building 5 Area, PR

Percent Solids:

1	File ID	DF	Analyzed	Ву	Prep Date	Prep Batch	Analytical Batch
Run #1	3M62832.D	1	07/07/16	LK	06/24/16	OP95059A	E3M2967

Run #2

Initial Volume Final Volume Run #1

Run #2

900 ml 1.0 ml

IXUII # E						
CAS No.	Compound	Result	RL	MDL	Units	Ç
91-20-3	Naphthalene	1.09	0.11	0.033	ug/l	
123-91-1	1,4-Dioxane	0.568	0.11	0.054	ug/l	
CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Lim	its	
4165-60-0	Nitrobenzene-d5	76%		24-1	25%	
321-60-8	2-Fluorobiphenyl	66%		19-1	27%	
1718-51-0	Terphenyl-d14	100%		10-1	19%	
	• •					



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

Report of Analysis

Page 1 of 1

Client Sample ID: MW-5 Lab Sample ID: JC22613-2

Matrix: Method:

AQ - Ground Water SW846-8015C (DAI)

Date Sampled: 06/17/16 Date Received: 06/21/16

Percent Solids: n/a

Project: BMSMC, Building 5 Area, PR

Run #1	File ID	DF	Analyzed	By	Prep Date	Prep Batch	Analytical Batch
Run #2	GH105635.D	1	06/23/16	XPL	n/a	n/a	GGH5333

Low Molecular Alcohol List

CAS No.	Compound	Result	RL	MDL	Units	Q
64-17-5 78-83-1 67-63-0 71-23-8 71-36-3 78-92-2 67-56-1	Ethanol Isobutyl Alcohol Isopropyl Alcohol n-Propyl Alcohol n-Butyl Alcohol sec-Butyl Alcohol	ND ND ND ND ND ND	100 100 100 100 100 100	55 36 68 43 87 66	ug/l ug/l ug/l ug/l ug/l ug/l	
CAS No.	Methanol Surrogate Recoveries Hexanol	ND Run# 1 102%	200 Run# 2	71 Lim 56-1	ug/I its 45%	



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

Page 1 of 3

E3E3716

SGS Accutest

Report of Analysis

SB

Client Sample ID: MW-13 Lab Sample ID: JC22613-3

File ID

Matrix: Method:

Project:

AQ - Ground Water SW846 8270D SW846 3510C

BMSMC, Building 5 Area, PR

07/02/16

Date Sampled: 06/17/16 Date Received: 06/21/16

Percent Solids: n/a

OP95059

Q

DF Analyzed **Analytical Batch** By Prep Date Prep Batch 3E84768.D

06/24/16

Run #1 Run #2

> Initial Volume Final Volume 910 ml

Run #1 Run #2

1.0 ml

1

ABN TCL Special List

CAS No.	Compound	Result	RL	MDL	Units
95-57-8	2-Chlorophenol	ND	5.5	0.90	ug/l
59-50-7	4-Chloro-3-methyl phenol	ND	5.5	0.98	ug/l
120-83-2	2,4-Dichlorophenol	ND	2.2	1.4	ug/l
105-67-9	2,4-Dimethylphenol	ND	5.5	2.7	ug/l
51-28-5	2,4-Dinitrophenol	ND	11	1.7	ug/l
534-52-1	4,6-Dinitro-o-cresol	ND	5.5	1.4	ug/l
95-48-7	2-Methylphenol	ND	2.2	0.98	ug/l
	3&4-Methylphenol	ND	2.2	0.97	ug/l
88-75-5	2-Nitrophenol	ND	5.5	1.1	ug/l
100-02-7	4-Nitrophenol	ND	11	1.3	ug/l
87-86-5	Pentachlorophenol	ND	5.5	1.5	ug/I
108-95-2	Phenol	ND	2.2	0.43	ug/l
58-90-2	2,3,4,6-Tetrachlorophenol	ND	5.5	1.6	ug/l
95-95-4	2,4,5-Trichlorophenol	ND	5.5	1.5	ug/l
88-06-2	2,4,6-Trichlorophenol	ND	5.5	1.0	ug/l
83-32-9	Acenaphthene	ND	1.1	0.21	ug/l
208-96-8	Acenaphthylene	ND	1.1	0.15	ug/l
98-86-2	Acetophenone	ND	2.2	0.23	ug/l
120-12-7	Anthracene	ND	1.1	0.23	ug/i
1912-24-9	Atrazine	ND	2.2	0.49	ug/l
100-52-7	Benzaldehyde	ND	5.5	0.32	ug/l
56-55-3	Benzo(a)anthracene	ND	1.1	0.22	ug/l
50-32-8	Benzo(a)pyrene	ND	1.1	0.23	ug/l
205-99-2	Benzo(b)fluoranthene	ND	1.1	0.23	ug/l
191-24-2	Benzo(g,h,i)perylene	ND	1.1	0.37	ug/l
207-08-9	Benzo(k)fluoranthene	ND	1.1	0.23	ug/l
101-55-3	4-Bromophenyl phenyl ether	ND	2.2	0.44	ug/l
85-68-7	Butyl benzyl phthalate	ND	2.2	0.50	ug/l
92-52-4	1,1'-Biphenyl	ND	1.1	0.23	ug/l
91-58-7	2-Chloronaphthalene	ND	2.2	0.26	ug/l
106-47-8	4-Chloroaniline	1.8	5.5	0.37	ug/l
86-74-8	Carbazole	ND	1.1	0.25	ug/I

fael Infante Méndez 10 = 1888

ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

N = Indicates presumptive evidence of a compound

J

Project:

Report of Analysis

Client Sample ID: MW-13 Lab Sample ID: JC22613-3

Matrix: Method:

AQ - Ground Water

SW846 8270D SW846 3510C BMSMC, Building 5 Area, PR

Date Sampled: 06/17/16 Date Received: 06/21/16

Percent Solids:

ABN TCL Special List

CAS No.	Compound	Result	RL	MDL	Units	Q
105-60-2	Caprolactam	ND	2.2	0.71	ug/l	
218-01-9	Chrysene	ND	1.1	0.19	ug/l	
111-91-1	bis (2-Chloroethoxy) methane	ND	2.2	0.31	ug/l	
111-44-4	bis(2-Chloroethyl)ether	ND	2.2	0.27	ug/l	
108-60-1	bis(2-Chloroisopropyl)ether	ND	2.2	0.44	ug/l	
7005-72-3	4-Chlorophenyl phenyl ether	ND	2.2	0.40	ug/i	
121-14-2	2,4-Dinitrotoluene	ND	1.1	0.61	ug/l	
606-20-2	2,6-Dinitrotoluene	ND	1.1	0.52	ug/l	
91-94-1	3,3'-Dichlorobenzidine	ND	2.2	0.56	ug/l	
53-70-3	Dibenzo(a,h)anthracene	ND	1.1	0.36	ug/l	
132-64-9	Dibenzofuran	ND	5.5	0.24	ug/l	
84-74-2	Di-n-butyl phthalate	ND	2.2	0.55	ug/l	
117-84-0	Di-n-octyl phthalate	ND	2.2	0.26	ug/l	
84-66-2	Diethyl phthalate	ND	2.2	0.29	ug/I	
131-11-3	Dimethyl phthalate	ND	2.2	0.24	ug/l	
117-81-7	bis(2-Ethylhexyl)phthalate	ND	2.2	1.8	ug/l	
206-44-0	Fluoranthene	0.71	1.1	0.19	ug/l	J
86-73-7	Fluorene	ND	1.1	0.19	ug/l	,
118-74-1	Hexachlorobenzene	ND	1.1	0.36	ug/l	
87-68-3	Hexachlorobutadiene	ND	1.1	0.54	ug/l	
77-47-4	Hexachlorocyclopentadiene	ND	11	3.1	ug/l	
67-72-1	Hexachloroethane	ND	2.2	0.43	ug/l	
193-39-5	Indeno(1,2,3-cd)pyrene	ND	1.1	0.36	ug/l	
78-59-1	Isophorone	ND	2.2	0.30	ug/l	
90-12-0	1-Methylnaphthalene	ND	1.1	0.29	ug/l	
91-57-6	2-Methylnaphthalene	ND	1.1	0.23	ug/l	
88-74-4	2-Nitroaniline	ND	5.5	0.30	ug/l	
99-09-2	3-Nitroaniline	ND	5.5	0.43	ug/l	
100-01-6	4-Nitroaniline	ND	5.5	0.48	ug/l	
98-95-3	Nitrobenzene	ND	2.2	0.71	ug/l	
621-64-7	N-Nitroso-di-n-propylamine	ND	2.2	0.53	ug/l	
86-30-6	N-Nitrosodiphenylamine	ND	5.5	0.24	ug/l	7
85-01-8	Phenanthrene	ND	1.1	0.19	ug/l	/:
129-00-0	Pyrene	0.50	1.1	0.24	ug/l	J (i
95-94-3	1,2,4,5-Tetrachlorobenzene	ND	2.2	0.41	ug/l	1
CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Lim	its	/
367-12-4	2-Fluorophenol	52%		14-8	8%	
4165-62-2	Phenol-d5	34%		10-1	10 %	



MDL = Method Detection Limit

RL = Reporting Limit

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

N = Indicates presumptive evidence of a compound

fuel Infante Méndez IC = 1888

E = Indicates value exceeds calibration range

Client Sample ID: MW-13 Lab Sample ID:

JC22613-3 AQ - Ground Water

Date Sampled: 06/17/16

Matrix: Method:

SW846 8270D SW846 3510C

Date Received: 06/21/16

Project:

BMSMC, Building 5 Area, PR

Percent Solids: n/a

ABN TCL Special List

CAS No.	Surrogate Recoveries	Run#1	Run# 2	Limits
118-79-6	2,4,6-Tribromophenol	58%		39-149%
4165-60-0	Nitrobenzene-d5	78%		32-128%
321-60-8	2-Fluorobiphenyl	77%		35-119%
1718-51-0	Terphenyl-d14	93%		10-126%



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

Report of Analysis

Page 1 of 1

Client Sample ID: Lab Sample ID:

MW-13 JC22613-3

Matrix:

AQ - Ground Water SW846 8270D BY SIM SW846 3510C Date Sampled: 06/17/16

Q

Date Received: 06/21/16

Method: Project:

BMSMC, Building 5 Area, PR

Percent Solids: n/a

File ID DF Ву Analyzed Prep Date Prep Batch **Analytical Batch** Run #1 3M62833.D 06/24/16 1 07/07/16 LK OP95059A E3M2967

DI

MOT

Run #2

Run #2

CASNo

Initial Volume Run #1 910 ml

Final Volume

Compound

1.0 ml

CAS No.	Compound	Result	RL	MDL	Units
91-20-3	Naphthalene	ND	0.11	0.032	ug/l
123-91-1	1,4-Dioxane	0.152	0.11	0.054	ug/l

CAS No.	Surrogate Recoveries	Run#1	Run# 2	Limits
4165-60-0	Nitrobenzene-d5	68%		24-125%
321-60-8	2-Fluorobiphenyl	63%		19-127%
1718-51-0	Terphenyl-d14	94%		10-119%



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

Report of Analysis

Page 1 of 1

Client Sample ID: MW-13 Lab Sample ID:

JC22613-3

Matrix: Method: AQ - Ground Water

Project:

SW846-8015C (DAI)

BMSMC, Building 5 Area, PR

Date Sampled: 06/17/16

Date Received: 06/21/16

Percent Solids: n/a

File ID DF Analyzed By Run #1 GH105636.D 1 06/23/16 XPL Run #2	Prep Date Prep Batch n/a n/a	Analytical Batch GGH5333
--	------------------------------	-----------------------------

Low Molecular Alcohol List

CAS No.	Compound	Result	RL	MDL	Units	Q
64-17-5	Ethanol	ND	100	55	ug/l	
78-83-1	Isobutyl Alcohol	ND	100	36	ug/l	
67-63-0	Isopropyl Alcohol	ND	100	68	ug/l	
71-23-8	n-Propyl Alcohol	ND	100	43	ug/l	
71-36-3	n-Butyl Alcohol	ND	100	87	ug/l	
78-92-2	sec-Butyl Alcohol	ND	100	66	ug/l	
67-56-1	Methanol	ND	200	71	ug/l	
CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Lim	its	
111-27-3	Hexanol	107%		56-1	45%	



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

N = Indicates presumptive evidence of a compound

Report of Analysis

Page 1 of 3

Client Sample ID: BR-4

Lab Sample ID: JC22613-4

Matrix: Method: AQ - Ground Water

Date Sampled: 06/17/16

Q

SW846 8270D SW846 3510C

DF

1

Date Received: 06/21/16

Project:

BMSMC, Building 5 Area, PR

Percent Solids: n/a

Run #1

File ID 3E84769.D

Analyzed By 07/02/16 SB Prep Date 06/24/16

Prep Batch OP95059

Analytical Batch E3E3716

Run #2

Initial Volume

1000 ml

Final Volume 1.0 ml

Run #1

Run #2

ABN TCL Special List

CAS No.	Compound	Result	RL	MDL	Unite
95-57-8	2-Chlorophenol	ND	5.0	0.82	ug/l
59-50-7	4-Chloro-3-methyl phenol	ND	5.0	0.89	ug/l
120-83-2	2,4-Dichlorophenol	ND	2.0	1.3	ug/i
105-67-9	2,4-Dimethylphenol	ND	5.0	2.4	ug/l
51-28-5	2,4-Dinitrophenol	ND	10	1.6	ug/l
534-52-1	4,6-Dinitro-o-cresol	ND	5.0	1.3	ug/l
95-48-7	2-Methylphenol	ND	2.0	0.89	ug/l
	3&4-Methylphenol	ND	2.0	0.88	ug/l
88-75-5	2-Nitrophenol	ND	5.0	0.96	ug/l
100-02-7	4-Nitrophenol	ND	10	1.2	ug/l
87-86-5	Pentachlorophenol	ND	5.0	1.4	ug/l
108-95-2	Phenol	ND	2.0	0.39	ug/l
58-90-2	2,3,4,6-Tetrachlorophenol	ND	5.0	1.5	ug/l
95-95-4	2,4,5-Trichlorophenol	ND	5.0	1.3	ug/l
88-06-2	2,4,6-Trichlorophenol	ND	5.0	0.92	ug/l
83-32-9	Acenaphthene	ND	1.0	0.19	ug/l
208-96-8	Acenaphthylene	ND	1.0	0.14	ug/l
98-86-2	Acetophenone	ND	2.0	0.21	ug/l
120-12-7	Anthracene	ND	1.0	0.21	ug/l
1912-24-9	Atrazine	ND	2.0	0.45	ug/i
100-52-7	Benzaldehyde	ND	5.0	0.29	ug/l
56-55-3	Benzo(a)anthracene	ND	1.0	0.20	ug/l
50-32-8	Benzo(a)pyrene	ND	1.0	0.21	ug/l
205-99-2	Benzo(b)fluoranthene	ND	1.0	0.21	ug/l
191-24-2	Benzo(g,h,i)perylene	ND	1.0	0.34	ug/l
207-08-9	Benzo(k)fluoranthene	ND	1.0	0.21	ug/l
101-55-3	4-Bromophenyl phenyl ether	ND	2.0	0.40	ug/l
85-68-7	Butyl benzyl phthalate	ND	2.0	0.46	ug/l
92-52-4	1,1'-Biphenyl	ND	1.0	0.21	ug/l
91-58-7	2-Chloronaphthalene	ND	2.0	0.24	ug/l
106-47-8	4-Chloroaniline	ND	5.0	0.34	ug/l
86-74-8	Carbazole	ND	1.0	0.23	ug/l
					_

iael Infante Méndez

ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

Client Sample ID: BR-4 Lab Sample ID: Matrix:

JC22613-4

AQ - Ground Water

Date Sampled: 06/17/16 Date Received: 06/21/16

Q

Method: Project:

SW846 8270D SW846 3510C BMSMC, Building 5 Area, PR Percent Solids: n/a

ABN TCL Special List

CAS No.	Compound	Result	RL	MDL	Units
105-60-2	Caprolactam	ND	2.0	0.65	ug/l
218-01-9	Chrysene	ND	1.0	0.18	ug/l
111-91-1	bis(2-Chloroethoxy)methane	ND	2.0	0.28	ug/l
111-44-4	bis (2-Chloroethyl) ether	ND	2.0	0.25	ug/f
108-60-1	bis (2-Chloroisopropyl) ether	ND	2.0	0.40	ug/l
7005-72-3	4-Chlorophenyl phenyl ether	ND	2.0	0.37	ug/l
121-14-2	2,4-Dinitrotoluene	ND	1.0	0.55	ug/l
606-20-2	2,6-Dinitrotoluene	ND	1.0	0.48	ug/l
91-94-1	3,3'-Dichlorobenzidine	ND	2.0	0.51	ug/l
53-70-3	Dibenzo(a,h)anthracene	ND	1.0	0.33	ug/l
132-64-9	Dibenzofuran	ND	5.0	0.22	ug/l
84-74-2	Di-n-butyl phthalate	ND	2.0	0.50	ug/l
117-84-0	Di-n-octyl phthalate	ND	2.0	0.23	ug/l
84-66-2	Diethyl phthalate	ND	2.0	0.26	ug/l
131-11-3	Dimethyl phthalate	ND	2.0	0.22	ug/l
117-81-7	bis(2-Ethylhexyl)phthalate	ND	2.0	1.7	ug/l
206-44-0	Fluoranthene	ND	1.0	0.17	ug/l
86-73-7	Fluorene	ND	1.0	0.17	ug/l
118-74-1	Hexachlorobenzene	ND	1.0	0.33	ug/l
87-68-3	Hexachlorobutadiene	ND	1.0	0.49	ug/l
77-47-4	Hexachlorocyclopentadiene	ND	10	2.8	ug/l
67-72-1	Hexachloroethane	ND	2.0	0.39	ug/l
193-39-5	Indeno(1,2,3-cd)pyrene	ND	1.0	0.33	ug/l
78-59-1	Isophorone	ND	2.0	0.28	ug/l
90-12-0	1-Methylnaphthalene	ND	1.0	0.26	ug/l
91-57-6	2-Methylnaphthalene	ND	1.0	0.21	ug/l
88-74-4	2-Nitroaniline	ND	5.0	0.28	ug/l
99-09-2	3-Nitroaniline	ND	5.0	0.39	ug/l
100-01-6	4-Nitroaniline	ND	5.0	0.44	ug/l
98-95-3	Nitrobenzene	ND	2.0	0.64	ug/l
621-64-7	N-Nitroso-di-n-propylamine	ND	2.0	0.48	ug/l
86-30-6	N-Nitrosodiphenylamine	ND	5.0	0.22	ug/l
85-01-8	Phenanthrene	ND	1.0	0.18	ug/l
129-00-0	Pyrene	ND	1.0	0.22	ug/l
95-94-3	1,2,4,5-Tetrachlorobenzene	ND	2.0	0.37	ug/l
CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Lim	its
367-12-4	2-Fluorophenol	48%		14-8	8%
4165-62-2	Phenol-d5	31%		10-1	10%



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

Client Sample ID: BR-4 Lab Sample ID:

JC22613-4

Matrix: Method:

Project:

AQ - Ground Water

SW846 8270D SW846 3510C

BMSMC, Building 5 Area, PR

Date Sampled: 06/17/16 Date Received: 06/21/16

Percent Solids: n/a

ABN TCL Special List

CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limits
118-79-6	2,4,6-Tribromophenol	58%		39-149%
4165-60-0	Nitrobenzene-d5	80%		32-128%
321-60-8	2-Fluorobiphenyl	77%		35-119%
1718-51-0	Terphenyl-d14	80%		10-126%



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

Report of Analysis

Page 1 of 1

Client Sample ID: BR-4

Lab Sample ID:

JC22613-4 Matrix:

File ID

3M62834.D

Method:

AQ - Ground Water

DF

1

SW846 8270D BY SIM SW846 3510C

Date Sampled:

Q

06/17/16 Date Received: 06/21/16

Percent Solids: n/a

Project:

BMSMC, Building 5 Area, PR

Analyzed By Prep Date Prep Batch Analytical Batch 07/07/16 LK 06/24/16 OP95059A E3M2967

10-119%

Run #1 Run #2

Initial Volume Final Volume Run #1 1000 ml 1.0 ml

Terphenyl-d14

Run #2

1718-51-0

CAS No. Compound Result RL MDL Units 91-20-3 Naphthalene ND 0.10 0.029 ug/l 123-91-1 1,4-Dioxane 1.49 0.10 0.049 ug/I CAS No. Surrogate Recoveries Run#1 Run# 2 Limits 4165-60-0 Nitrobenzene-d5 66% 24-125% 19-127% 321-60-8 2-Fluorobiphenyl 69%

78%



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

Page 1 of 1

SGS Accutest

Report of Analysis

Client Sample ID: BR-4

Lab Sample ID: JC22613-4

Matrix: Method: AQ - Ground Water SW846-8015C (DAI)

Project:

BMSMC, Building 5 Area, PR

Date Sampled: 06/17/16

Date Received: 06/21/16

Percent Solids: n/a

Run #1	File ID	DF	Analyzed	By	Prep Date	Prep Batch	Analytical Batch
Run #2	GH105637.D	1	06/23/16	XPL	n/a	n/a	GGH5333

Low Molecular Alcohol List

CAS No.	Compound	Result	RL	MDL	Units	Q
64-17-5 78-83-1 67-63-0 71-23-8 71-36-3 78-92-2	Ethanol Isobutyl Alcohol Isopropyl Alcohol n-Propyl Alcohol n-Butyl Alcohol sec-Butyl Alcohol	ND ND ND ND ND	100 100 100 100 100	55 36 68 43 87	ug/l ug/l ug/l ug/l ug/l ug/l	
67-56-1	Methanol	ND	200	71	ug/l	
CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Lim	its	
111-27-3	Hexanol	107%		56-1	45%	





MDL = Method Detection Limit

RL = Reporting Limit

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

E = Indicates value exceeds calibration range

Report of Analysis

Page 1 of 3

Client Sample ID: MW-16

Lab Sample ID: JC22613-5

Matrix: Method: AQ - Ground Water

SW846 8270D SW846 3510C

Date Sampled: 06/17/16 Date Received: 06/21/16

Percent Solids: n/a

Project: BMSMC, Building 5 Area, PR

File ID DF Analyzed By Prep Date Run #1 07/02/16 3E84770.D 1 SĎ 06/24/16

Prep Batch **OP95059**

Q

Analytical Batch E3E3716

Run #2

Initial Volume Final Volume

Run #1

990 ml 1.0 ml

Run #2

ABN TCL Special List

CAS No.	Compound	Result	RL	MDL	Unit
95-57-8	2-Chlorophenol	ND	5.1	0.83	ug/l
59-50-7	4-Chloro-3-methyl phenol	ND	5.1	0.90	ug/l
120-83-2	2,4-Dichlorophenol	ND	2.0	1.3	ug/l
105-67-9	2,4-Dimethylphenol	ND	5.1	2.5	ug/l
51-28-5	2,4-Dinitrophenol	ND	10	1.6	ug/I
534-52-1	4,6-Dinitro-o-cresol	ND	5.1	1.3	ug/l
95-48-7	2-Methylphenol	ND	2.0	0.90	ug/l
	3&4-Methylphenol	ND	2.0	0.89	ug/l
88-75-5	2-Nitrophenol	ND	5.1	0.97	ug/l
100-02-7	4-Nitrophenol	ND	10	1.2	ug/l
87-86-5	Pentachlorophenol	ND	5.1	1.4	ug/l
108-95-2	Phenol	ND	2.0	0.40	ug/l
58-90-2	2,3,4,6-Tetrachlorophenol	ND	5.1	1.5	ug/i
95-95-4	2,4,5-Trichlorophenol	ND	5.1	1.3	ug/l
88-06-2	2,4,6-Trichlorophenol	ND	5.1	0.93	ug/l
83-32-9	Acenaphthene	ND	1.0	0.19	ug/l
208-96-8	Acenaphthylene	ND	1.0	0.14	ug/l
98-86-2	Acetophenone	ND	2.0	0.21	ug/l
120-12-7	Anthracene	ND	1.0	0.21	ug/l
1912-24-9	Atrazine	ND	2.0	0.45	ug/l
100-52-7	Benzaldehyde	ND	5.1	0.29	ug/l
56-55-3	Benzo(a)anthracene	ND	1.0	0.21	ug/l
50-32-8	Benzo(a) pyrene	ND	1.0	0.22	ug/l
205-99-2	Benzo(b)fluoranthene	ND	1.0	0.21	ug/l
191-24-2	Benzo(g,h,i)perylene	ND	1.0	0.34	ug/l
207-08-9	Benzo(k)fluoranthene	ND	1.0	0.21	ug/l
101-55-3	4-Bromophenyl phenyl ether	ND	2.0	0.41	ug/l
85-68-7	Butyl benzyl phthalate	ND	2.0	0.46	ug/l
92-52-4	1,1'-Biphenyl	ND	1.0	0.21	ug/l
91-58-7	2-Chloronaphthalene	ND	2.0	0.24	ug/l
106-47-8	4-Chloroaniline	ND	5.1	0.34	ug/l
86-74-8	Carbazole	ND	1.0	0.23	ug/I
					3



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

Report of Analysis

Client Sample ID:	MW-16
Lab Sample ID:	JC22613-5

Matrix: AQ - Ground Water

Method: SW846 8270D SW846 3510C Project:

BMSMC, Building 5 Area, PR

Date Sampled: 06/17/16 Date Received: 06/21/16

Percent Solids: n/a

Q

ABN	TCL	Special	List
-----	-----	----------------	------

CAS No.	Compound	Result	RL	MDL	Units
105-60-2	Caprolactam	ND	2.0	0.66	ug/l
218-01-9	Chrysene	ND	1.0	0.18	ug/l
111-91-1	bis(2-Chloroethoxy)methane	ND	2.0	0.28	ug/l
111-44-4	bis (2-Chloroethyl) ether	ND	2.0	0.25	ug/l
108-60-1	bis(2-Chloroisopropyl)ether	ND	2.0	0.41	ug/l
7005-72-3	4-Chlorophenyl phenyl ether	ND	2.0	0.37	ug/l
121-14-2	2,4-Dinitrotoluene	ND	1.0	0.56	ug/l
606-20-2	2,6-Dinitrotoluene	ND	1.0	0.48	ug/l
91-94-1	3,3'-Dichlorobenzidine	ND	2.0	0.51	ug/l
53-70-3	Dibenzo(a,h)anthracene	ND	1.0	0.33	ug/l
132-64-9	Dibenzofuran	ND	5.1	0.22	ug/l
84-74-2	Di-n-butyl phthalate	ND	2.0	0.50	ug/l
117-84-0	Di-n-octyl phthalate	ND	2.0	0.24	ug/l
84-66-2	Diethyl phthalate	ND	2.0	0.26	ug/l
131-11-3	Dimethyl phthalate	ND	2.0	0.22	ug/l
117-81-7	bis(2-Ethylhexyl)phthalate	ND	2.0	1.7	ug/l
206-44-0	Fluoranthene	ND	1.0	0.17	ug/l
86-73-7	Fluorene	ND	1.0	0.17	ug/l
118-74-1	Hexachlorobenzene	ND	1.0	0.33	ug/l
87-68-3	Hexachlorobutadiene	ND	1.0	0.50	ug/l
77-47-4	Hexachlorocyclopentadiene	ND	10	2.8	ug/l
67-72-1	Hexachloroethane	ND	2.0	0.39	ug/l
193-39-5	Indeno(1,2,3-cd)pyrene	ND	1.0	0.34	ug/l
78-59-1	Isophorone	ND	2.0	0.28	ug/l
90-12-0	1-Methylnaphthalene	ND	1.0	0.27	ug/l
91-57-6	2-Methylnaphthalene	ND	1.0	0.21	ug/l
88-74-4	2-Nitroaniline	ND	5.1	0.28	ug/i
99-09-2	3-Nitroaniline	ND	5.1	0.39	ug/l
100-01-6	4-Nitroaniline	ND	5.1	0.44	ug/l
98-95-3	Nitrobenzene	ND	2.0	0.65	ug/l
621-64-7	N-Nitroso-di-n-propylamine	ND	2.0	0.49	ug/l
86-30-6	N-Nitrosodiphenylamine	ND	5.1	0.22	ug/l
85-01-8	Phenanthrene	ND	1.0	0.18	ug/l
129-00-0	Pyrene	ND	1.0	0.22	ug/l
95-94-3	1,2,4,5-Tetrachlorobenzene	ND	2.0	0.37	ug/l
CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limi	ts
367-12-4	2-Fluorophenol	46%		14-8	3%
4165-62-2	Phenol-d5	29%		10-1	10%

fact infante Méndez

ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

Client Sample ID: Lab Sample ID:

MW-16 JC22613-5

AQ - Ground Water SW846 8270D SW846 3510C Date Sampled: Date Received: 06/21/16 Percent Solids: n/a

06/17/16

Method: Project:

Matrix:

BMSMC, Building 5 Area, PR

ABN TCL Special List

CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limits
118-79-6	2,4,6-Tribromophenol	60%		39-149%
4165-60-0	Nitrobenzene-d5	79%		32-128%
321-60-8	2-Fluorobiphenyl	82%		35-119%
1718-51-0	Terphenyl-d14	97%		10-126%



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

Report of Analysis

Page 1 of 1

Client Sample ID:	MW-16
Lab Sample ID:	JC22613-

Matrix: Method:

Project:

AQ - Ground Water

SW846 8270D BY SIM SW846 3510C

BMSMC, Building 5 Area, PR

Date Sampled: 06/17/16

Q

Date Received: 06/21/16

Percent Solids: n/a

File ID DF Ву **Analytical Batch** Analyzed Prep Date Prep Batch Run #1 3M62835.D 07/07/16 LK 06/24/16 OP95059A E3M2967

Run #2

Initial Volume Final Volume Run #1 990 ml 1.0 ml

Run #2

CAS No.	Compound	Result	RL	MDL	Units	- (
91-20-3 123-91-1	Naphthalene 1,4-Dioxane	ND 0.334	0.10 0.10	0.030 0.049	ug/l ug/l	
CAS No.	Surrogate Recoveries	Run#1	Run# 2	Limi	its	
4165-60-0	Nitrobenzene-d5	68%		24-1	25%	
321-60-8	2-Fluorobiphenyl	68%		19-1	27%	
1718-51-0	Terphenyl-d14	94%		10-1		



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

Project:

Report of Analysis

Page 1 of 1

Client Sample ID: MW-16 Lab Sample ID:

JC22613-5

Matrix: Method:

AQ - Ground Water

SW846-8015C (DAI)

BMSMC, Building 5 Area, PR

Date Sampled:

06/17/16 Date Received: 06/21/16

Percent Solids: n/a

File ID DF **Analytical Batch** Analyzed By Prep Date Prep Batch Run #1 GH105638.D XPL 06/23/16 GGH5333 n/a n/a Run #2

Low Molecular Alcohol List

CAS No.	Compound	Result	RL	MDL	Units	Q
64-17-5	Ethanol	ND	100	55	ug/l	
78-83-1	Isobutyl Alcohol	ND	100	36	ug/l	
67-63-0	Isopropyl Alcohol	ND	100	68	ug/l	
71-23-8	n-Propyl Alcohol	ND	100	43	ug/l	
71-36-3	n-Butyl Alcohol	ND	100	87	ug/l	
78-92-2	sec-Butyl Alcohol	ND	100	66	ug/l	
67-56-1	Methanol	ND	200	71	ug/l	
CAS No.	Surrogate Recoveries	Run#1	Run# 2	Lim	its	
111-27-3	Hexanol	110%		56-1	45%	



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

N = Indicates presumptive evidence of a compound

Report of Analysis

Page 1 of 3

Client Sample ID: MW-9 Lab Sample ID: JC22613-6

Matrix: Method: AQ - Ground Water

SW846 8270D SW846 3510C

Date Sampled: 06/17/16 Date Received: 06/21/16

Percent Solids: n/a

Q

Project: BMSMC, Building 5 Area, PR

File ID DF Analyzed Ву **Analytical Batch** Prep Date Prep Batch Run #1 3E84771.D 1 07/02/16 SB 06/24/16 OP95059 E3E3716

Run #2

Initial Volume Final Volume 1000 ml 1.0 ml

Run #1 Run #2

ABN TCL Special List

CAS No.	Compound	Result	RL	MDL	Units
95-57-8	2-Chlorophenol	ND	5.0	0.82	ug/l
59-50-7	4-Chloro-3-methyl phenol	ND	5.0	0.89	ug/l
120-83-2	2,4-Dichlorophenol	ND	2.0	1.3	ug/l
105-67-9	2,4-Dimethylphenol	ND	5.0	2.4	ug/l
51-28-5	2,4-Dinitrophenol	ND	10	1.6	ug/l
534-52-1	4,6-Dinitro-o-cresol	ND	5.0	1.3	ug/l
95-48-7	2-Methylphenol	ND	2.0	0.89	ug/l
	3&4-Methylphenol	ND	2.0	0.88	ug/l
88-75-5	2-Nitrophenol	ND	5.0	0.96	ug/l
100-02-7	4-Nitrophenol	ND	10	1.2	ug/l
87-86-5	Pentachlorophenol	ND	5.0	1.4	ug/l
108-95-2	Phenol	ND	2.0	0.39	ug/l
58-90-2	2,3,4,6-Tetrachlorophenol	ND	5.0	1.5	ug/l
95-95-4	2,4,5-Trichlorophenol	ND	5.0	1.3	ug/l
88-06-2	2,4,6-Trichlorophenol	ND	5.0	0.92	ug/l
83-32-9	Acenaphthene	ND	1.0	0.19	ug/l
208-96-8	Acenaphthylene	ND	1.0	0.14	ug/l
98-86-2	Acetophenone	ND	2.0	0.21	ug/l
120-12-7	Anthracene	ND	1.0	0.21	ug/i
1912-24-9	Atrazine	ND	2.0	0.45	ug/l
100-52-7	Benzaldehyde	ND	5.0	0.29	ug/l
56-55-3	Benzo(a)anthracene	ND	1.0	0.20	ug/l
50-32-8	Benzo(a) pyrene	ND	1.0	0.21	ug/l
205-99-2	Benzo(b)fluoranthene	ND	1.0	0.21	ug/l
191-24-2	Benzo(g,h,i)perylene	ND	1.0	0.34	ug/l
207-08-9	Benzo(k)fluoranthene	ND	1.0	0.21	ug/l
101-55-3	4-Bromophenyl phenyl ether	ND	2.0	0.40	ug/l
85-68-7	Butyl benzyl phthalate	ND	2.0	0.46	ug/l
92-52-4	1,1'-Biphenyl	ND	1.0	0.21	ug/l
91-58-7	2-Chloronaphthalene	ND	2.0	0.24	ug/l
106-47-8	4-Chloroaniline	ND	5.0	0.34	ug/l
86-74-8	Carbazole	ND	1.0	0.23	ug/i

fuel Infan Méndez IC # 1881

ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank



Report of Analysis

Client Sample ID: MW-9 Lab Sample ID: JC22613-6

Matrix: AQ - Ground Water

Method: SW846 8270D SW846 3510C Project:

BMSMC, Building 5 Area, PR

Date Sampled: 06/17/16 Date Received: 06/21/16

Percent Solids: n/a

Q

ABN TCL Special List

CAS No.	Compound	Result	RL	MDL	Units
105-60-2	Caprolactam	ND	2.0	0.65	ug/l
218-01-9	Chrysene	ND	1.0	0.18	ug/l
111-91-1	bis(2-Chloroethoxy)methane	ND	2.0	0.28	ug/l
111-44-4	bis(2-Chloroethyl)ether	ND	2.0	0.25	ug/l
108-60-1	bis (2-Chloroisopropyl) ether	ND	2.0	0.40	ug/l
7005-72-3	4-Chlorophenyl phenyl ether	ND	2.0	0.37	ug/i
121-14-2	2,4-Dinitrotoluene	ND	1.0	0.55	ug/l
606-20-2	2,6-Dinitrotoluene	ND	1.0	0.48	ug/i
91-94-1	3,3'-Dichlorobenzidine	ND	2.0	0.51	ug/l
53-70-3	Dibenzo(a,h)anthracene	ND	1.0	0.33	ug/l
132-64-9	Dibenzofuran	ND	5.0	0.22	ug/l
84-74-2	Di-n-butyl phthalate	ND	2.0	0.50	ug/l
117-84-0	Di-n-octyl phthalate	ND	2.0	0.23	ug/l
84-66-2	Diethyl phthalate	ND	2.0	0.26	ug/l
131-11-3	Dimethyl phthalate	ND	2.0	0.22	ug/l
117-81-7	bis(2-Ethylhexyl)phthalate	ND	2.0	1.7	ug/l
206-44-0	Fluoranthene	ND	1.0	0.17	ug/l
86-73-7	Fluorene	ND	1.0	0.17	ug/l
118-74-1	Hexachlorobenzene	ND	1.0	0.33	ug/l
87-68-3	Hexachlorobutadiene	ND	1.0	0.49	ug/l
77-47-4	Hexachlorocyclopentadiene	ND	10	2.8	ug/l
67-72-1	Hexachloroethane	ND	2.0	0.39	ug/l
193-39-5	Indeno(1,2,3-cd)pyrene	ND	1.0	0.33	ug/l
78-59-1	Isophorone	ND	2.0	0.28	ug/l
90-12-0	1-Methylnaphthalene	ND	1.0	0.26	ug/l
91-57-6	2-Methylnaphthalene	ND	1.0	0.21	ug/l
88-74-4	2-Nitroaniline	ND	5.0	0.28	ug/l
99-09-2	3-Nitroaniline	ND	5.0	0.39	ug/l
100-01-6	4-Nitroaniline	ND	5.0	0.44	ug/l
98-95-3	Nitrobenzene	ND	2.0	0.64	ug/l
621-64-7	N-Nitroso-di-n-propylamine	ND	2.0	0.48	ug/l
86-30-6	N-Nitrosodiphenylamine	ND	5.0	0.22	ug/i
85-01-8	Phenanthrene	ND	1.0	0.18	ug/l
129-00-0	Pyrene	ND	1.0	0.22	ug/I
95-94-3	1,2,4,5-Tetrachlorobenzene	ND	2.0	0.37	ug/l
CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Lim	its
367-12-4	2-Fluorophenol	52%		14-8	8%
4165-62-2	Phenol-d5	33%		10-1	10%



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

N = Indicates presumptive evidence of a compound

Method:

Project:

Client Sample ID: MW-9 Lab Sample ID: Matrix:

JC22613-6

AQ - Ground Water

SW846 8270D SW846 3510C BMSMC, Building 5 Area, PR

Date Sampled: 06/17/16 Date Received: 06/21/16

Percent Solids: n/a

ABN TCL Special List

CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limits
118-79-6	2,4,6-Tribromophenol	61%		39-149%
4165-60-0	Nitrobenzene-d5	83%		32-128%
321-60-8	2-Fluorobiphenyl	81%		35-119%
1718-51-0	Terphenyl-d14	93%		10-126%



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

Report of Analysis

Page 1 of 1

Client Sample ID: Lab Sample ID:

MW-9 JC22613-6

Matrix: Method: AQ - Ground Water

SW846 8270D BY SIM SW846 3510C

Date Sampled: Date Received: 06/21/16

06/17/16

BMSMC, Building 5 Area, PR

Percent Solids: n/a

Q

1	File ID	DF	Analyzed	Bv	Prep Date	Prep Batch	Analytical Batch
1	_			-,	1. op 2000	Trop Decem	1 may aver Daver
Run #1	3M62836.D	1	07/07/16	T K	06/24/16	OP95059A	E3M2967
12.000	DIVIOLOGO.L	-	01/01/10	LIC	00/24/10	OI 22022V	POINTPOOL
Run #2							

(un #2

Project:

Run #1 Run #2

Initial Volume Final Volume 1.0 ml

1000 ml

CAS No.	Compound	Result	RL	MDL	Units

91-20-3 Naphthalene ND 0.10 0.029ug/l 123-91-1 1,4-Dioxane 0.398 0.10 0.049ug/l CAS No. Surrogate Recoveries Run#1 Run# 2 Limits

4165-60-0 Nitrobenzene-d5 71% 24-125% 321-60-8 2-Fluorobiphenyl 72% 19-127% 1718-51-0 Terphenyl-d14 92% 10-119%



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

Report of Analysis

Page 1 of 1

Client Sample ID: Lab Sample ID:

MW-9 JC22613-6

Matrix: Method: AQ - Ground Water SW846-8015C (DAI)

DF

1

Date Sampled:

06/17/16 Date Received: 06/21/16

Percent Solids: n/a

Project:

BMSMC, Building 5 Area, PR

Run #1 Run #2

File ID

GH105639.D

Analyzed 06/23/16

By XPL

Prep Batch Prep Date

Analytical Batch GGH5333

Low Molecular Alcohol List

CAS No.	Compound	Result	RL	MDL	Units	Q
64-17-5	Ethanol	ND	100	55	ug/l	
78-83-1	Isobutyl Alcohol	ND	100	36	ug/l	
67-63-0	Isopropyl Alcohol	ND	100	68	ug/l	
71-23-8	n-Propyl Alcohol	ND	100	43	ug/l	
71-36-3	n-Butyl Alcohol	ND	100	87	ug/l	
78-92-2	sec-Butyl Alcohol	ND	100	66	ug/l	
67-56-1	Methanol	ND	200	71	ug/l	
CAS No.	Surrogate Recoveries	Run#1	Run# 2	Lim	its	
111-27-3	Hexanol	101%		56-1	45%	



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

Report of Analysis

Page 1 of 3

Client Sample ID: Lab Sample ID:

MW-23S JC22613-7

Matrix: Method: AQ - Ground Water

SW846 8270D SW846 3510C

Date Sampled: Date Received: 06/21/16

06/20/16

Percent Solids: n/a

Project:

BMSMC, Building 5 Area, PR

Run #1

File ID 3E84772.D DF Analyzed 1 07/02/16

By Prep Date SB 06/24/16

Prep Batch OP95059

Q

Analytical Batch E3E3716

Run #2

Initial Volume Final Volume

Run #1

840 ml

1.0 ml

Run #2

ABN TCL Special List

CAS No.	Compound	Result	RL	MDL	Unit
95-57-8	2-Chlorophenol	ND	6.0	0.98	ug/l
59-50-7	4-Chloro-3-methyl phenol	ND	6.0	1.1	ug/l
120-83-2	2,4-Dichlorophenol	ND	2.4	1.5	ug/l
105-67-9	2,4-Dimethylphenol	ND	6.0	2.9	ug/l
51-28-5	2,4-Dinitrophenol	ND	12	1.8	ug/l
534-52-1	4,6-Dinitro-o-cresol	ND	6.0	1.5	ug/l
95-48-7	2-Methylphenol	ND	2.4	1.1	ug/l
	3&4-Methylphenol	ND	2.4	1.0	ug/l
88-75-5	2-Nitrophenol	ND	6.0	1.1	ug/l
100-02-7	4-Nitrophenol	ND	12	1.4	ug/l
87-86-5	Pentachlorophenol	ND	6.0	1.6	ug/l
108-95-2	Phenol	ND	2.4	0.47	ug/l
58-90-2	2,3,4,6-Tetrachlorophenol	ND	6.0	1.7	ug/l
95-95-4	2,4,5-Trichlorophenol	ND	6.0	1.6	ug/l
88-06-2	2,4,6-Trichlorophenol	ND	6.0	1,1	ug/l
83-32-9	Acenaphthene	ND	1.2	0.23	ug/l
208-96-8	Acenaphthylene	ND	1.2	0.16	ug/l
98-86-2	Acetophenone	ND	2.4	0.25	ug/l
120-12-7	Anthracene	ND	1.2	0.25	ug/l
1912-24-9	Atrazine	ND	2.4	0.53	ug/l
100-52-7	Benzaldehyde	ND	6.0	0.34	ug/l
56-55-3	Benzo(a)anthracene	ND	1.2	0.24	ug/l
50-32-8	Benzo(a) pyrene	ND	1.2	0.25	ug/l
205-99-2	Benzo(b) fluoranthene	ND	1.2	0.24	ug/l
191-24-2	Benzo(g,h,i)perylene	ND	1.2	0.41	ug/l
207-08-9	Benzo(k)fluoranthene	ND	1.2	0.24	ug/l
101-55-3	4-Bromophenyl phenyl ether	ND	2.4	0.48	ug/l
85-68-7	Butyl benzyl phthalate	ND	2.4	0.54	ug/l
92-52-4	1,1'-Biphenyi	ND	1.2	0.25	ug/l
91-58-7	2-Chloronaphthalene	ND	2.4	0.28	ug/l
106-47-8	4-Chloroaniline	ND	6.0	0.40	ug/l
86-74-8	Carbazole	ND	1.2	0.27	ug/l
					_



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank



Client Sample ID: Lab Sample ID:

MW-23S

Matrix:

JC22613-7 AQ - Ground Water

Method:

SW846 8270D SW846 3510C

Project:

BMSMC, Building 5 Area, PR

Date Sampled: 06/20/16

Date Received: 06/21/16

Percent Solids: n/a



ABN TCL Special List

ABN TCL	Special List					
CAS No.	Compound	Result	RL	MDL	Units	Q
105-60-2	Caprolactam	ND	2.4	0.77	ug/l	
218-01-9	Chrysene	ND	1.2	0.21	ug/l	
111-91-1	bis (2-Chloroethoxy) methane	ND	2.4	0.33	ug/l	
111-44-4	bis (2-Chloroethyl) ether	ND	2.4	0.30	ug/l	
108-60-1	bis (2-Chloroisopropyl)ether	ND	2.4	0.48	ug/l	
7005-72-3	4-Chlorophenyl phenyl ether	ND	2.4	0.44	ug/l	
121-14-2	2,4-Dinitrotoluene	ND	1.2	0.66	ug/l	
606-20-2	2,6-Dinitrotoluene	ND	1.2	0.57	ug/l	
91-94-1	3,3'-Dichlorobenzidine	ND	2.4	0.60	ug/l	
53-70-3	Dibenzo(a,h)anthracene	ND	1.2	0.39	ug/l	
132-64-9	Dibenzofuran	ND	6.0	0.26	ug/l	
84-74-2	Di-n-butyl phthalate	ND	2.4	0.59	ug/l	
117-84-0	Di-n-octyl phthalate	ND	2.4	0.28	ug/l	
84-66-2	Diethyl phthalate	ND	2.4	0.31	ug/l	
131-11-3	Dimethyl phthalate	ND	2.4	0.26	ug/l	
117-81-7	bis(2-Ethylhexyl)phthalate	ND	2.4	2.0	ug/l	
206-44-0	Fluoranthene	ND	1.2	0.20	ug/l	
86-73-7	Fluorene	ND	1.2	0.20	ug/l	
118-74-1	Hexachlorobenzene	ND	1.2	0.39	ug/l	
87-68-3	Hexachlorobutadiene	ND	1.2	0.59	ug/l	
77-47-4	Hexachlorocyclopentadiene	ND	12	3.3	ug/l	
67-72-1	Hexachloroethane	ND	2.4	0.46	ug/l	
193-39-5	Indeno(1,2,3-cd)pyrene	ND	1.2	0.40	ug/l	
78-59-1	Isophorone	ND	2.4	0.33	ug/i	
90-12-0	1-Methylnaphthalene	ND	1.2	0.31	ug/l	
91-57-6	2-Methylnaphthalene	ND	1.2	0.25	ug/l	
88-74-4	2-Nitroaniline	ND	6.0	0.33	ug/l	
99-09-2	3-Nitroaniline	ND	6.0	0.46	ug/l	
100-01-6	4-Nitroaniline	ND	6.0	0.52	ug/l	
98-95-3	Nitrobenzene	ND	2.4	0.76	ug/l	
621-64-7	N-Nitroso-di-n-propylamine	ND	2.4	0.57	ug/l	
86-30-6	N-Nitrosodiphenylamine	ND	6.0	0.26	ug/i	200
85-01-8	Phenanthrene	ND	1.2	0.21	ug/l	SOCIALO DE
129-00-0	Pyrene	ND	1.2	0.26	ug/l	all land
95-94-3	1,2,4,5-Tetrachlorobenzene	ND	2.4	0.44	ug/l	
					Ū	fael Infante
CAS No.	Surrogate Recoveries	Run#1	Run# 2	Lim	its	Méndez
						IC = 1888
367-12-4	2-Fluorophenol	57%		14-8	8%	1/1
4165-62-2	Phenol-d5	36%		10-1	10%	MICO LICENCIA

ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank



Client Sample ID: MW-23S Lab Sample ID: Matrix:

JC22613-7

AQ - Ground Water

Date Sampled: 06/20/16 Date Received: 06/21/16 Percent Solids: n/a

Method: Project:

SW846 8270D SW846 3510C BMSMC, Building 5 Area, PR

ABN TCL Special List

CAS No.	Surrogate Recoveries	Run#1	Run# 2	Limits
118-79-6	2,4,6-Tribromophenol	59%		39-149%
4165-60-0	Nitrobenzene-d5	81%		32-128%
321-60-8	2-Fluorobiphenyl	84%		35-119%
1718-51-0	Terphenyl-d14	98%		10-126%



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

Report of Analysis

By

LK

Page 1 of 1

Client Sample ID: MW-23S Lab Sample ID: JC22613-7

Matrix: Method: AQ - Ground Water

SW846 8270D BY SIM SW846 3510C

Date Sampled: 06/20/16 Date Received: 06/21/16

Percent Solids: n/a

Project: BMSMC, Building 5 Area, PR

File ID DF Analyzed Run #1 3M62837.D 07/07/16

Prep Date 06/24/16

Prep Batch OP95059A

Analytical Batch E3M2967

Initial Volume Final Volume

Run #1 Run #2

Run #2

840 ml

1.0 ml

CAS No.	Compound	Result	RL	MDL Units		Q
91-20-3 123-91-1	Naphthalene 1,4-Dioxane	ND ND	0.12 0.12	0.035 0.058	ug/l ug/l	
				Limits		
CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Lim	its	
CAS No. 4165-60-0	Surrogate Recoveries Nitrobenzene-d5	Run# 1 71%	Run# 2		i ts 25%	
	-		Run# 2	24-1		



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

Report of Analysis

Page 1 of 1

Client Sample ID: MW-23S Lab Sample ID: JC22613-7

Matrix: Method:

Project:

AQ - Ground Water

SW846-8015C (DAI)

BMSMC, Building 5 Area, PR

Date Sampled: Date Received:

06/20/16 06/21/16

Percent Solids: n/a

Run #1	File ID	DF	Analyzed	By	Prep Date	Prep Batch	Analytical Batch
	GH105642.D	1	06/23/16	XPL	n/a	n/a	GGH5333
Run #2							

Low Molecular Alcohol List

CAS No.	Compound	Result	RL	MDL	Units	Q
64-17-5	Ethanol	ND	100	55	ug/l	
78-83-1	Isobutyl Alcohol	ND	100	36	ug/l	
67-63-0	Isopropyl Alcohol	ND	100	68	ug/l	
71-23-8	n-Propyl Alcohol	ND	100	43	ug/l	
71-36-3	n-Butyl Alcohol	ND	100	87	ug/l	
78-92-2	sec-Butyl Alcohol	ND	100	66	ug/l	
67-56-1	Methanol	ND	200	71	ug/i	
CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limits		
111-27-3	Hexanol	102%	102%		56-145%	



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

N = Indicates presumptive evidence of a compound

SGS Accutest

Report of Analysis

By

DS

Page 1 of 1

Client Sample ID: MW-23S Lab Sample ID:

JC22613-7

Date Sampled: 06/20/16

Matrix:

AQ - Ground Water

DF

1

Date Received: 06/21/16

Method:

SW846 8081B SW846 3510C

Percent Solids: n/a

Project:

BMSMC, Building 5 Area, PR

Run #1 Run #2 File ID 6G36640.D Analyzed 06/27/16

Prep Date 06/27/16

Prep Batch OP95079

Analytical Batch G6G1047

Initial Volume **Final Volume**

1000 ml

10.0 ml

Run #1 Run #2

Pesticide TCL List

CAS No.	Compound	Result	RL	MDL	Units
309-00-2	Aldrin	ND	0.010	0.0060	ug/l
319-84-6	alpha-BHC	ND	0.010	0.0060	ug/l
319-85-7	beta-BHC	ND	0.010	0.0057	ug/l
319-86-8	delta-BHC	ND	0.010	0.0046	ug/l
58-89-9	gamma-BHC (Lindane)	ND	0.010	0.0028	ug/l
5103-71-9	alpha-Chlordane	ND	0.010	0.0046	ug/l
5103-74-2	gamma-Chlordane	ND	0.010	0.0046	ug/l
60-57-1	Dieldrin	ND	0.010	0.0036	ug/l
72-54-8	4,4'-DDD	ND	0.010	0.0038	ug/l
72-55-9	4,4'-DDE	ND	0.010	0.0062	ug/l
50-29-3	4,4'-DDT	ND	0.010	0.0050	ug/l
72-20-8	Endrin	ND	0.010	0.0050	ug/l
1031-07-8	Endosulfan sulfate	ND	0.010	0.0053	ug/l
7421-93-4	Endrin aldehyde	ND	0.010	0.0051	ug/l
53494-70-5	Endrin ketone	ND	0.010	0.0051	ug/l
959-98-8	Endosulfan-I	ND	0.010	0.0050	ug/l
33213-65-9	Endosulfan-II	ND	0.010	0.0043	ug/l
76-44-8	Heptachlor	ND	0.010	0.0038	ug/l
1024-57-3	Heptachlor epoxide	ND	0.010	0.0065	ug/l
72-43-5	Methoxychlor	ND	0.020	0.0057	ug/l
8001-35-2	Toxaphene	ND	0.25	0.18	ug/l
CAS No.	Surrogate Recoveries	Run#1	Run# 2	Limi	ts
877-09-8	Tetrachloro-m-xylene	76%		26-13	32%
877-09-8	Tetrachloro-m-xylene	77%		26-13	32%
2051-24-3	Decachlorobiphenyl	83%		10-1	18%
2051-24-3	Decachlorobiphenyl	91%		10-1	18%



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

SGS Accutest

Report of Analysis

Page 1 of 3

Client Sample ID: MW-21S

Lab Sample ID: JC22613-8

File ID

Matrix: Method: AQ - Ground Water

SW846 8270D SW846 3510C

Date Sampled: Date Received:

Q

06/20/16 06/21/16

Percent Solids: n/a

Project:

BMSMC, Building 5 Area, PR

Run #1 Run #2

DF 3E84773.D 1

Analyzed 07/02/16

By Prep Date SB 06/24/16

Prep Batch OP95059

Analytical Batch

E3E3716

Initial Volume 910 ml

Final Volume 1.0 ml

Run #1

Run #2

ABN TCL Special List

CAS No.	Compound	Result	RL	MDL	Unita
95-57-8	2-Chlorophenol	ND	5.5	0.90	ug/l
59-50-7	4-Chloro-3-methyl phenol	ND	5.5	0.98	ug/l
120-83-2	2,4-Dichlorophenol	ND	2.2	1.4	ug/l
105-67-9	2,4-Dimethylphenol	ND	5.5	2.7	ug/l
51-28-5	2,4-Dinitrophenol	ND	11	1.7	ug/l
534-52-1	4,6-Dinitro-o-cresol	ND	5.5	1.4	ug/l
95-48-7	2-Methylphenol	ND	2.2	0.98	ug/l
	3&4-Methylphenol	ND	2.2	0.97	ug/l
88-75-5	2-Nitrophenol	ND	5.5	1.1	ug/l
100-02-7	4-Nitrophenol	ND	11	1.3	ug/l
87-86-5	Pentachlorophenol	ND	5.5	1.5	ug/l
108-95-2	Phenol	ND	2.2	0.43	ug/l
58-90-2	2,3,4,6-Tetrachlorophenol	ND	5.5	1.6	ug/l
95-95-4	2,4,5-Trichlorophenol	ND	5.5	1.5	ug/l
88-06-2	2,4,6-Trichlorophenol	ND	5.5	1.0	ug/l
83-32-9	Acenaphthene	ND	1.1	0.21	ug/I
208-96-8	Acenaphthylene	ND	1.1	0.15	ug/l
98-86-2	Acetophenone	ND	2.2	0.23	ug/I
120-12-7	Anthracene	ND	1.1	0.23	ug/l
1912-24-9	Atrazine	ND	2.2	0.49	ug/l
100-52-7	Benzaldehyde	ND	5.5	0.32	ug/l
56-55-3	Benzo(a)anthracene	ND	1.1	0.22	ug/l
50-32-8	Benzo(a)pyrene	ND	1.1	0.23	ug/l
205-99-2	Benzo(b)fluoranthene	ND	1.1	0.23	ug/l
191-24-2	Benzo(g,h,i)perylene	ND	1.1	0.37	ug/l
207-08-9	Benzo(k)fluoranthene	ND	1.1	0.23	ug/l
101-55-3	4-Bromophenyl phenyl ether	ND	2.2	0.44	ug/l
85-68-7	Butyl benzyl phthalate	ND	2.2	0.50	ug/l
92-52-4	1,1'-Biphenyl	ND	1.1	0.23	ug/l
91-58-7	2-Chloronaphthalene	ND	2.2	0.26	ug/I
106-47-8	4-Chloroaniline	ND	5.5	0.37	ug/l
86-74-8	Carbazole	ND	1.1	0.25	ug/l



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

Report of Analysis

Page 2 of 3

Client Sample ID: MW-21S

Lab Sample ID: JC22613-8

Matrix: Method:

Project:

AQ - Ground Water

SW846 8270D SW846 3510C BMSMC, Building 5 Area, PR Date Sampled: 06/20/16 Date Received: 06/21/16

Percent Solids: n/a

Q

ABN TCL Special List

CAS No.	Compound	Result	RL,	MDL	Units
105-60-2	Caprolactam	ND	2.2	0.71	ug/l
218-01-9	Chrysene	ND	1.1	0.19	ug/l
111-91-1	bis(2-Chloroethoxy)methane	ND	2.2	0.31	ug/l
111-44-4	bis(2-Chloroethyl)ether	ND	2.2	0.27	ug/l
108-60-1	bis (2-Chloroisopropyl) ether	ND	2.2	0.44	ug/l
7005-72-3	4-Chlorophenyl phenyl ether	ND	2.2	0.40	ug/l
121-14-2	2,4-Dinitrotoluene	ND	1.1	0.61	ug/l
606-20-2	2,6-Dinitrotoluene	ND	1.1	0.52	ug/l
91-94-1	3,3'-Dichlorobenzidine	ND	2.2	0.56	ug/l
53-70-3	Dibenzo(a,h)anthracene	ND	1.1	0.36	ug/l
132-64-9	Dibenzofuran	ND	5.5	0.24	ug/l
84-74-2	Di-n-butyl phthalate	ND	2.2	0.55	ug/l
117-84-0	Di-n-octyl phthalate	ND	2.2	0.26	ug/l
84-66-2	Diethyl phthalate	ND	2.2	0.29	ug/l
131-11-3	Dimethyl phthalate	ND	2.2	0.24	ug/l
117-81-7	bis (2-Ethylhexyl) phthalate	ND	2.2	1.8	ug/l
206-44-0	Fluoranthene	ND	1.1	0.19	ug/l
86-73-7	Fluorene	ND	1.1	0.19	ug/l
118-74-1	Hexachlorobenzene	ND	1.1	0.36	ug/l
87-68-3	Hexachlorobutadiene	ND	1.1	0.54	ug/l
77-47-4	Hexachlorocyclopentadiene	ND	11	3.1	ug/l
67-72-1	Hexachloroethane	ND	2.2	0.43	ug/l
193-39-5	Indeno(1,2,3-cd)pyrene	ND	1.1	0.36	ug/l
78-59-1	Isophorone	ND	2.2	0.30	ug/I
90-12-0	1-Methylnaphthalene	ND	1.1	0.29	ug/l
91-57-6	2-Methylnaphthalene	ND	1.1	0.23	ug/l
88-74-4	2-Nitroaniline	ND	5.5	0.30	ug/l
99-09-2	3-Nitroaniline	ND	5.5	0.43	ug/l
100-01-6	4-Nitroaniline	ND	5.5	0.48	ug/l
98-95-3	Nitrobenzene	ND	2.2	0.71	ug/l
621-64-7	N-Nitroso-di-n-propylamine	ND	2.2	0.53	ug/l
86-30-6	N-Nitrosodiphenylamine	ND	5.5	0.24	ug/l
85-01-8	Phenanthrene	ND	1.1	0.19	ug/l
129-00-0	Pyrene	ND	1.1	0.24	ug/l
95-94-3	1,2,4,5-Tetrachlorobenzene	ND	2.2	0.41	ug/I
CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limi	ts
367-12-4	2-Fluorophenol	53%		14-88	3%
4165-62-2	Phenol-d5	34%		10-1	10%



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

MW-21S

Client Sample ID: Lab Sample ID: JC22613-8

Matrix:

AQ - Ground Water

Date Sampled: Date Received: 06/21/16

06/20/16

Percent Solids: n/a

Method: Project:

SW846 8270D SW846 3510C BMSMC, Building 5 Area, PR

Report of Analysis

ABN TCL Special List

CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limits
118-79-6	2,4,6-Tribromophenol	57%		39-149%
4165-60-0	Nitrobenzene-d5	79%		32-128%
321-60-8	2-Fluorobiphenyl	80%		35-119%
1718-51-0	Terphenyl-d14	87%		10-126%





MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank



1718-51-0

Terphenyl-d14

Report of Analysis

Page 1 of 1

Client Sam Lab Sampl Matrix: Method: Project:	E ID: JC2261 AQ - G SW846	3-8 round Wate 8270D BY		3510C		Date		5/20/16 5/21/16 /a
Run #1 Run #2	File ID 3M62838.D	DF 1	Analyzed 07/07/16	By LK	Prep D 06/24/1		Prep Batch OP95059A	Analytical Batch E3M2967
Run #1 Run #2	Initial Volume 910 ml	Final Vo	lume					
CAS No.	Compound		Result	RL	MDL	Units	Q	
91-20-3 123-91-1	Naphthalene 1,4-Dioxane		ND 0.646	0.11 0.11	0.032 0.054	ug/l ug/l		
CAS No.	Surrogate Rec	overies	Run# 1	Run# 2	Lim	its		
4165-60-0 321-60-8	Nitrobenzene-d 2-Fluorobipher		73% 73%			.25% .27%		

90%



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

10-119%

B = Indicates analyte found in associated method blank

SGS Accutest

Report of Analysis

Page 1 of 1

Client Sample ID: MW-21S

Lab Sample ID: JC22613-8 Matrix:

AQ - Ground Water SW846-8015C (DAI) Date Sampled: Date Received:

06/20/16 06/21/16

Method: Project:

BMSMC, Building 5 Area, PR

Percent Solids: n/a

Run #1	File ID	DF	Analyzed	By	Prep Date	Prep Batch	Analytical Batch
	GH105643.D	1	06/23/16	XPL	n/a	n/a	GGH5333
Run #2							

Low Molecular Alcohol List

CAS No.	Compound	Result	RL	MDŁ	Units	Q
64-17-5	Ethanol	ND	100	55	ug/l	
78-83-1	Isobutyl Alcohol	ND	100	36	ug/l	
67-63-0	Isopropyl Alcohol	ND	100	68	ug/l	
71-23-8	n-Propyl Alcohol	ND	100	43	ug/l	
71-36-3	n-Butyl Alcohol	ND	100	87	ug/l	
78-92-2	sec-Butyl Alcohol	ND	100	66	ug/I	
67-56-1	Methanol	ND	200	71	ug/l	
CAS No.	Surrogate Recoveries	Run#1	Run# 2	Lim	its	
111-27-3	Hexanol	99%		56-1	45%	



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

SGS Accutest

Report of Analysis

Page 1 of 1

Client Sample ID: MW-21S

Lab Sample ID: JC22613-8

Matrix: Method:

Project:

AQ - Ground Water

SW846 8081B SW846 3510C

Date Sampled: 06/20/16 Date Received: 06/21/16

Q

Percent Solids: n/a

BMSMC, Building 5 Area, PR

File ID DF Analyzed By Prep Date Run #1 6G36643.D 06/28/16 DS 06/27/16 1

Prep Batch **Analytical Batch** OP95079 G6G1047

Run #2

Initial Volume Final Volume

Run #1

980 ml

10.0 ml

Run #2

Pesticide TCL List

CAS No.	Compound	Result	RL	MDL	Units
309-00-2	Aldrin	ND	0.010	0.0062	ug/l
319-84-6	alpha-BHC	ND	0.010	0.0061	ug/l
319-85-7	beta-BHC	ND	0.010	0.0058	ug/I
319-86-8	delta-BHC	ND	0.010	0.0047	ug/l
58-89-9	gamma-BHC (Lindane)	ND	0.010	0.0028	ug/l
5103-71-9	alpha-Chlordane	ND	0.010	0.0047	ug/l
5103-74-2	gamma-Chlordane	ND	0.010	0.0047	ug/l
60-57-1	Dieldrin	ND	0.010	0.0037	ug/l
72-54-8	4,4'-DDD	ND	0.010	0.0039	ug/l
72-55-9	4,4'-DDE	ND	0.010	0.0063	ug/l
50-29-3	4,4'-DDT	ND	0.010	0.0051	ug/l
72-20-8	Endrin	ND	0.010	0.0051	ug/l
1031-07-8	Endosulfan sulfate	ND	0.010	0.0054	ug/l
7421-93-4	Endrin aldehyde	ND	0.010	0.0052	ug/l
53494-70-5	Endrin ketone	ND	0.010	0.0052	ug/l
959-98-8	Endosulfan-I	ND	0.010	0.0051	ug/l
33213-65-9	Endosulfan-II	ND	0.010	0.0044	ug/l
76-44-8	Heptachlor	ND	0.010	0.0039	ug/l
1024-57-3	Heptachlor epoxide	ND	0.010	0.0067	ug/l
72-43-5	Methoxychlor	ND	0.020	0.0058	ug/l
8001-35-2	Toxaphene	ND	0.26	0.19	ug/l
CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limi	ts
877-09-8	Tetrachloro-m-xylene	91%		26-13	32%
877-09-8	Tetrachloro-m-xylene	90%		26-13	32%
2051-24-3	Decachlorobiphenyl	88%		10-11	18%
2051-24-3	Decachlorobiphenyl	96%		10-11	18%



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

SGS ACC	UTEST -			SGS Ac Louise 130 9-020u	natest - De	ijton , NJ 10681 12-329-34	0				900 14	-	 P59			-		ter Carrie		/ 元		613
nderson Mulhalland Ass	Project Name:	BMS	2	(在)	1	500	und	IA AC	ter				8				250		10			Matrix Code
	Direct	ر ۱۰۰		3		010				550	2700	2	1808	हर्					9		- 11	W - Ground W WW - Water W - Surface W
700 Westchester	cq.		-34	Conye	y Name	e (II allber		n Repar	t loj	<u> </u>	\dashv ²	0151	8	24							-	SCI - Soil SL- Skinge
inchase NY	Hum	1040	PR	Shreet A								9	3	ŽΊ				- }			Ш	SED-Sedmer OI - DII UQ - Direc Leg
Cerry Taylor	Promot if			Server 4	37908							-	3	39							- 1	ASR - Arr SOL - Other S
14-251-0100	Chard Purchase	Order if		Cay			Shee	io .		Zφ		1.4	J.	2.4						ì		WP - Wps FB-Flow Star
estor Rivera	Proad Harase			Alterdica	,			_					빙	뙲					9	ı		Equament (P.S Pinee Bin 120-Tito Bien
estor Nivera	+	T	Culturan	1		т	660			Holms.	غ ⊢	7	5	× 4					ĝ.	1	L	140-170 Dan
				Sample]		_ [8	10000 10000	y à	8 5	SVOC	Ę	Posticides)A	1	Γ	
Field ID / Point of Collection	SECHIONAL I	aliali.	fire and and	-	Tables 1		₽ §	1 D	Because of	보호	_	+=	0	_		-				\dashv	+	LAB USE OF
MW 3	+	6/16/16	1545	NR	GW	_	3	-	2	H	X	×	\vdash	쏫	\dashv	H	\dashv	-			-	€9 V46
MW-5 MW-13	-	617116	1242	NR.	GW GW		3	+	2	 	×	∤	\vdash	×		\vdash				\dashv	-	V 10 1
RR - 4	-	117/11	1505	NA	6W		3	H	1		- \(\)	₹ V		分	\vdash					-	7	
MW-16		6/17/16	1602	NR	6w	_	9		2	Ш	X	×		X								
MW 9		6/17/16	1934	NR	GW		3		2		×	X		Y							1	
MW-235		6/20/16	1235		CW		3		3	Ш	<u> </u> \	X,	X	X						_	4	
Mw 215	-	12014	1547	NA	6N	6	3		3_		-X	X	X	x			\dashv	\dashv		\dashv	+	
							\blacksquare		H		+	F					\Box			7	7	
Turnerund Erre Baumen days	SCORES N	9E 29E	-0-785	মূ-ন		Data	Destructi	data (rá	ormano.	Щ	関本	700	2 1973		\$13E}	Com	mards /	Space	i Instruc	ions 25	741	A. ·
1		Apparent Pile Date			Commerc	M "A" (L) M TIT (L)				YASP CA	Angury A		A	11 -6	9 (Cep.	4	1-1	teth	VI N	tal	thaler
Dat 18 Statemen Coys Day RUSH 1NITIAL A	SESSMENT	Am/B		E		Level 3rd			ᇊ	inte Fore	ne .		_			_				32	_	
1 2 Day RUSH LABEL V	ERIFICATION	T)X			Commerc	ad Tinoon	-		<u>_</u>													
1 Day RUSH					etad "A" =	Recurs On	у, Солг	الشاجة	B' - Re	eulin + OC	C Survey	y										
& Plant TIA man products VIA Lambon	- 8	emple Custody n	ust be done		niger and	in Nove sa	Symme mphes	ry + Per clustege	tul Ran poste	sae aalon, k	ncluding	courter		y.	-0	y is w	The K	upon	receip	in the	\$2	ratory
1100			-				2	7	=x					Date: To	2.1-1	6	2	484	\geq	7		_
		Services By:	<u> </u>				-	and by			. 2.20			n= 11	-		4	et Ber				Term on a
		13					-	_	_	_	■ 4900		-		-	_			-0			dela

. . .

JC22613: Chain of Custody Page 1 of 2

EXECUTIVE NARRATIVE

SDG No:

JC22613

Laboratory:

Accutest, New Jersey

Analysis:

SW846-8270D

Number of Samples:

8

Location:

BMSMC, Building 5 Area

Humacao, PR

SUMMARY: Eight (8) samples were analyzed for the ABN TCL list following method SW846-8270D; Naphthalene and 1,4-Dioxane were also analyzed by SW846-8270D using the selective ion monitoring (SIM) technique. The sample results were assessed according to USEPA data validation guidance documents in the following order of precedence: EPA Hazardous Waste Support Section, SOP HW-35A, July 2015 –Revision 0. Semivolatile Data Validation. The QC criteria and data validation actions listed on the data review worksheets are from the primary guidance document, unless otherwise noted.

Results are valid and can be used for decision making purposes.

Critical issues:

None

Major:

None

Minor:

None

Critical findings:

None

Major findings:

None

Minor findings:

 Initial and continuing calibration verifications meet the method and guidance document required performance criteria except in the cases described in the Data Review Worhseet. Analytes not meeting the continuing calibration verification method performance criteria and validation guidance document performance criteria qualified as estimated (J) or (UJ) in affected samples.

Analytes not meeting the continuing calibration verification method performance criteria but were within the validation guidance document performance criteria were not qualified. .

No closing calibration verification included in data package. No action taken, professional judgment.

QC samples from other jobs were not validated.

COMMENTS:

Results are valid and can be used for decision making purposes.

Reviewers Name:

Rafael Infante

Chemist License 1888

Signature:

Date:

July 21, 1016

SAMPLE ORGANIC DATA SAMPLE SUMMARY

Sample ID: JC22613-1

Sample location: BMSMC Building 5 Area

Sampling date: 6/16/2016 Matrix: Groundwater

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
2-Chlorophenol	5.6	ug/l	1	-	U	Yes
4-Chloro-3-methyl phenol	5.6	ug/l	1	-	U	Yes
2,4-Dichlorophenol	2.2	ug/l	1	•	U	Yes
2,4-Dimethylphenol	5.6	ug/l	1	-	U	Yes
2,4-Dinitrophenol	11	ug/l	1	-	U	Yes
4,6-Dinitro-o-cresol	5.6	ug/l	1	-	U	Yes
2-Methylphenol	2.2	ug/l	1	-	U	Yes
3&4-Methylphenol	2.2	ug/l	1	-	U	Yes
2-Nitrophenol	5.6	ug/l	1	-	U	Yes
4-Nitrophenol	11	ug/l	1	-	U	Yes
Pentachlorophenol	5.6	ug/l	1	-	U	Yes
Phenol	2.2	ug/l	1	-	U	Yes
2,3,4,6-Tetrachlorophenol	5.6	ug/l	1	_	UJ	Yes
2,4,5-Trichlorophenol	5.6	ug/l	1	-	U	Yes
2,4,6-Trichlorophenol	5.6	ug/l	1	-	U	Yes
Acenaphthene	0.50	ug/l	1	j	UJ	Yes
Acenaphthylene	1.1	ug/l	1	-	U	Yes
Acetophenone	2.2	ug/l	1	-	U	Yes
Anthracene	1.1	ug/l	1	-	U	Yes
Atrazine	2.2	ug/l	1	-	U	Yes
Benzaldehyde	5.6	ug/l	1	-	U	Yes
Benzo(a)anthracene	1.1	ug/l	1	-	U	Yes
Benzo(a)pyrene	1.1	ug/l	1	-	U	Yes
Benzo(b)fluoranthene	1.1	ug/l	1	-	U	Yes
Benzo(g,h,i)perylene	1.1	ug/l	1	-	U	Yes
Benzo(k)fluoranthene	1.1	ug/l	1	-	U	Yes
4-Bromophenyl phenyl ether	2.2	ug/l	1	-	U	Yes
Butyl benzyl phthalate	2.2	ug/l	1	-	U	Yes
1,1'-Biphenyl	1.1	ug/l	1	-	U	Yes
2-Chloronaphthalene	2.2	ug/l	1	-	U	Yes
4-Chloroaniline	5.6	ug/l	1	-	U	Yes
Carbazole	1.1	ug/l	1	-	U	Yes
Caprolactam	2.2	ug/l	1	-	U	Yes
Chrysene	1.1	ug/l	1	-	U	Yes
bis(2-Chloroethoxy)methane	2.2	ug/l	1	-	U	Yes
bis (2-Chloroethyl) ether	2.2	ug/l	1	-	U	Yes

bis(2-Chloroisopropyl)ether	2.2	ug/l	1		U	Yes
4-Chlorophenyl phenyl ether	2.2	ug/l	1	-	U	Yes
2,4-Dinitrotoluene	1.1	ug/l	1	-	U	Yes
2,6-Dinitrotoluene	1.1	ug/l	1	-	U	Yes
3,3'-Dichlorobenzidine	2.2	ug/l	1	-	U	Yes
Dibenzo(a,h)anthracene	1.1	ug/l	1	-	U	Yes
Dibenzofuran	5.6	ug/l	1		U	Yes
Di-n-butyl phthalate	2.2	ug/l	1		U	Yes
Di-n-octyl phthalate	2.2	ug/l	1	-	U	Yes
Diethyl phthalate	2.2	ug/l	1		U	Yes
Dimethyl phthalate	2.2	ug/l	1	-	U	Yes
bis(2-Ethylhexyl)phthalate	2.2	ug/l	1	-	U	Yes
Fluoranthene	1.1	ug/l	1		U	Yes
Fluorene	1.1	ug/l	1	-	-	Yes
Hexachlorobenzene	1.1	ug/l	1	-	U	Yes
Hexachlorobutadiene	1.1	ug/l	1	-	U	Yes
Hexachlorocyclopentadiene	11	ug/l	1	_	U	Yes
Hexachloroethane	2.2	ug/l	1	-	U	Yes
Indeno(1,2,3-cd)pyrene	1.1	ug/l	1		U	Yes
Isophorone	2.2	ug/l	1		U	Yes
1-Methylnaphthalene	35.9	ug/l	1	-	U	Yes
2-Methylnaphthalene	33.3	ug/l	1	-	Ų	Yes
2-Nitroaniline	5.6	ug/l	1	-	U	Yes
3-Nitroaniline	5.6	ug/l	1	-	U	Yes
4-Nitroaniline	5.6	ug/l	1	-	U	Yes
Nitrobenzene	2.2	ug/l	1	-	U	Yes
N-Nitroso-di-n-propylamine	2.2	ug/l	1	-	U	Yes
Nitrosodiphenylamine	5.6	ug/l	1	-	U	Yes
Phenanthrene	1.1	ug/l	1	-	U	Yes
Pyrene	1.1	ug/l	1	-	U	Yes
1,2,4,5-Tetrachlorobenzene	2.2	ug/l	1	-	U	Yes
METHOD:	8270D (SI	M)				
Naphthalene	2.19	ug/i	1	- 2	_	Yes
1,4-Dioxane	0.11	ug/l	1	-	U	Yes

2 B

Sample location: BMSMC Building 5 Area

Sampling date: 6/17/2016 Matrix: Groundwater

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
2-Chlorophenol	5.6	ug/l	1	-	U	Yes
4-Chloro-3-methyl phenol	5.6	ug/l	1	_	Ü	Yes
2,4-Dichlorophenol	2.2	ug/l	1	_	Ū	Yes
2,4-Dimethylphenol	5.6	ug/l	1	_	Ü	Yes
2,4-Dinitrophenol	11	ug/l	1	_	Ū	Yes
4,6-Dinitro-o-cresol	5.6	ug/l	1	_	Ü	Yes
2-Methylphenol	2.2	ug/l	1	_	U	Yes
3&4-Methylphenol	2.2	ug/i	1	_	Ü	Yes
2-Nitrophenol	5.6	ug/l	1	_	U	Yes
4-Nitrophenol	11	ug/l	1	-	U	Yes
Pentachlorophenol	5.6	ug/l	1	-	U	Yes
Phenol	2.2	ug/l	1	-	U	Yes
2,3,4,6-Tetrachlorophenol	5.6	ug/l	1	-	U	Yes
2,4,5-Trichlorophenol	5.6	ug/l	1	-	U	Yes
2,4,6-Trichlorophenol	5.6	ug/l	1	-	U	Yes
Acenaphthene	1.10	ug/l	1	-	U	Yes
Acenaphthylene	1.1	ug/l	1	-	U	Yes
Acetophenone	2.2	ug/l	1	-	U	Yes
Anthracene	1.1	ug/l	1	-	U	Yes
Atrazine	2.2	ug/l	1	-	U	Yes
Benzaldehyde	5.6	ug/l	1	-	U	Yes
Benzo(a)anthracene	1.1	ug/l	1	-	U	Yes
Benzo(a)pyrene	1.1	ug/l	1	-	U	Yes
Benzo(b)fluoranthene	1.1	ug/l	1	-	U	Yes
Benzo(g,h,i)perylene	1.1	ug/l	1	-	U	Yes
Benzo(k)fluoranthene	1.1	ug/l	1	-	U	Yes
4-Bromophenyl phenyl ether	2.2	ug/l	1	-	U	Yes
Butyl benzyl phthalate	2.2	ug/l	1	-	U	Yes
1,1'-Biphenyl	1.1	ug/i	1	~	U	Yes
2-Chloronaphthalene	2.2	ug/l	1	~	U	Yes
4-Chloroaniline	5.6	ug/l	1	-	U	Yes
Carbazole	1.1	ug/l	1	-	U	Yes
Caprolactam	2.2	ug/l	1	-	U	Yes
Chrysene	1.1	ug/l	1	-	U	Yes
bis(2-Chloroethoxy)methane	2.2	ug/l	1	-	U	Yes
bis(2-Chloroethyl)ether	2.2	ug/l	1	-	U	Yes
bis (2-Chlorois opropyl) ether	2.2	ug/l	1	-	U	Yes
4-Chlorophenyl phenyl ether	2.2	ug/l	1	-	U	Yes

2,4-Dinitrotoluene	1.1	ug/l	1	2	U	Yes
2,6-Dinitrotoluene	1.1	ug/l	1	+.	U	Yes
3,3'-Dichlorobenzidine	2.2	ug/l	1	-	U	Yes
Dibenzo(a,h)anthracene	1.1	ug/l	1	-	U	Yes
Dibenzofuran	5.6	ug/l	1	53	U	Yes
Di-n-butyl phthalate	2.2	ug/l	1	22	U	Yes
Di-n-octyl phthalate	2.2	ug/l	1	-	U	Yes
Diethyl phthalate	2.2	ug/l	1	-	U	Yes
Dimethyl phthalate	2.2	ug/l	1	-	U	Yes
bis(2-Ethylhexyl)phthalate	2.2	ug/l	1	100	U	Yes
Fluoranthene	1.1	ug/l	1	-	U	Yes
Fluorene	1.1	ug/l	1	2	U	Yes
Hexachlorobenzene	1.1	ug/l	1	1.5	U	Yes
Hexachlorobutadiene	1.1	ug/l	1	2	U	Yes
Hexachlorocyclopentadiene	11	ug/l	1	-	נט	Yes
Hexachloroethane	2.2	ug/l	1	÷:	U	Yes
Indeno(1,2,3-cd)pyrene	1.1	ug/l	1	-	U	Yes
Isophorone	2.2	ug/l	1	-	U	Yes
1-Methylnaphthalene	1.1	ug/l	1	-	U	Yes
2-Methylnaphthalene	1.1	ug/l	1	2.4	U	Yes
2-Nitroaniline	5.6	ug/l	1	-	U	Yes
3-Nitroaniline	5.6	ug/l	1	V 21	U	Yes
4-Nitroaniline	5.6	ug/l	1	-	U	Yes
Nitrobenzene	2.2	ug/l	1	-	U	Yes
N-Nitroso-di-n-propylamine	2.2	ug/l	1	_	UJ	Yes
Nitrosodiphenylamine	5.6	ug/l	1	-	U	Yes
Phenanthrene	1.1	ug/l	1	-	U	Yes
Pyrene	1.1	ug/l	Į	-	U	Yes
1,2,4,5-Tetrachlorobenzene	2.2	ug/l	1	-	U	Yes
AACTUOD	02700 /61	A 41				
METHOD:	827UD (SI	IVI)				
Naphthalene	1.09	ug/l	1	_	_	Yes
1,4-Dioxane	0.568	ug/l	1		-	Yes

Sample location: BMSMC Building 5 Area

Sampling date: 6/17/2016 Matrix: Groundwater

METHOD: 82/0D								
Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable		
2-Chlorophenol	5.5	ug/l	1	-	U	Yes		
4-Chloro-3-methyl phenol	5.5	ug/l	1	-	U	Yes		
2,4-Dichlorophenol	2.2	ug/l	1	-	U	Yes		
2,4-Dimethylphenol	5.5	ug/l	1	-	U	Yes		
2,4-Dinitrophenol	11	ug/l	1	-	U	Yes		
4,6-Dinitro-o-cresol	5.5	ug/l	1	-	U	Yes		
2-Methylphenol	2.2	ug/l	1	-	U	Yes		
3&4-Methylphenol	2.2	ug/l	1	-	U	Yes		
2-Nitrophenol	5.5	ug/l	1	-	U	Yes		
4-Nitrophenol	11	ug/l	1	-	U	Yes		
Pentachlorophenol	5.5	ug/l	1	-	U	Yes		
Phenol	2.2	ug/l	1	-	U	Yes		
2,3,4,6-Tetrachlorophenol	5.5	ug/l	1	-	U	Yes		
2,4,5-Trichlorophenol	5.5	ug/l	1	-	U	Yes		
2,4,6-Trichlorophenol	5.5	ug/l	1	-	U	Yes		
Acenaphthene	1.1	ug/l	1	-	U	Yes		
Acenaphthylene	1.1	ug/l	1	-	U	Yes		
Acetophenone	2.2	ug/l	1	-	U	Yes		
Anthracene	1.1	ug/l	1	-	U	Yes		
Atrazine	2.2	ug/l	1	-	U	Yes		
Benzaldehyde	5.5	ug/l	1	-	U	Yes		
Benzo(a)anthracene	1.1	ug/l	1	-	U	Yes		
Benzo(a)pyrene	1.1	ug/l	1	-	U	Yes		
Benzo(b)fluoranthene	1.1	ug/l	1	-	U	Yes		
Benzo(g,h,i)perylene	1.1	ug/l	1	-	U	Yes		
Benzo(k)fluoranthene	1.1	ug/l	1	-	U	Yes		
4-Bromophenyl phenyl ether	2.2	ug/l	1	-	U	Yes		
Butyl benzyl phthalate	2.2	ug/l	1	-	U	Yes		
1,1'-Biphenyl	1.1	ug/l	1	-	U	Yes		
2-Chloronaphthalene	2.2	ug/l	1	-	U	Yes		
4-Chloroaniline	5.5	ug/l	1	-	U	Yes		
Carbazole	1.1	ug/l	1	-	U	Yes		
Caprolactam	2.2	ug/l	1	-	U	Yes		
Chrysene	1.1	ug/l	1	-	U	Yes		
bis(2-Chloroethoxy)methane	2.2	ug/l	1	-	U	Yes		
bis(2-Chloroethyl)ether	2.2	ug/l	1	-	U	Yes		
bis (2-Chlorois opropyl) ether	2.2	ug/i	1	-	U	Yes		
4-Chlorophenyl phenyl ether	2.2	ug/l	1	-	U	Yes		

2,4-Dinitrotoluene	1.1	ug/l	1	-	U	Yes
2,6-Dinitrotoluene	1.1	ug/l	1	-	U	Yes
3,3'-Dichlorobenzidine	2.0	ug/l	1	_	U	Yes
Dibenzo(a,h)anthracene	1.1	ug/l	1	-	U	Yes
Dibenzofuran	5.5		1	-	U	
		ug/l	1	2.53	-	Yes
Di-n-butyl phthalate	2.2	ug/l	_	_	Ü	Yes
Di-n-octyl phthalate	2.2	ug/l	1	-	U	Yes
Diethyl phthalate	2.2	ug/l	1	-	U	Yes
Dimethyl phthalate	2.2	ug/l	1	-	U	Yes
bis(2-Ethylhexyl)phthalate	2.2	ug/l	1	-	U	Yes
Fluoranthene	0.71	ug/l	1	1	נט	Yes
Fluorene	1.1	ug/i	1	-	U	Yes
Hexachlorobenzene	1.1	ug/l	1		U	Yes
Hexachlorobutadiene	1.1	ug/l	1	720	U	Yes
Hexachlorocyclopentadiene	10	ug/l	1	-	UJ	Yes
Hexachloroethane	2.2	ug/l	1	-	U	Yes
Indeno(1,2,3-cd)pyrene	1.1	ug/l	1	_	U	Yes
Isophorone	2.0	ug/l	1	-	U	Yes
1-Methylnaphthalene	1.1	ug/l	1	_	U	Yes
2-Methylnaphthalene	1.1	ug/l	1	-	U	Yes
2-Nitroaniline	5.5	ug/l	1	-	U	Yes
3-Nitroaniline	5.5	ug/l	1		U	Yes
4-Nitroaniline	5.5	ug/l	1	-	U	Yes
Nitrobenzene	2.2	ug/l	1	-	U	Yes
N-Nitroso-di-n-propylamine	2.2	ug/l	1		Ū	Yes
Nitrosodiphenylamine	5.5	ug/l	_ 1	-	Ü	Yes
Phenanthrene	1.1	ug/l	1		Ü	Yes
Pyrene	1.1	ug/l	1	200	U	Yes
1,2,4,5-Tetrachlorobenzene	2.2	ug/l	1		U	Yes
1,2,7,3-1611 aci 1101 ODE1 126116	۷.۷	u6/1	1	3.5	U	162

METHOD: 8270D (SIM)

Naphthalene	0.11	ug/l	1	20	U	Yes
1,4-Dioxane	0.152	ug/l	1	-	-	Yes

Sample location: BMSMC Building 5 Area

Sampling date: 6/17/2016 Matrix: Groundwater

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
2-Chlorophenol	5.0	ug/i	1	-	U	Yes
4-Chloro-3-methyl phenol	5.0	ug/l	1	-	U	Yes
2,4-Dichlorophenol	2.0	ug/l	1	-	U	Yes
2,4-Dimethylphenol	5.0	ug/l	1	-	U	Yes
2,4-Dinitrophenol	10	ug/l	1	-	U	Yes
4,6-Dinitro-o-cresol	5.0	ug/l	1	-	U	Yes
2-Methylphenol	2.0	ug/l	1	-	U	Yes
3&4-Methylphenol	2.0	ug/l	1	-	U	Yes
2-Nitrophenol	5.0	ug/l	1	-	U	Yes
4-Nitrophenol	10	ug/l	1	-	U	Yes
Pentachlorophenol	5.0	ug/l	1	-	U	Yes
Phenol	2.0	ug/l	1	-	U	Yes
2,3,4,6-Tetrachlorophenol	5.0	ug/l	1	-	U	Yes
2,4,5-Trichlorophenol	5.0	ug/l	1	-	U	Yes
2,4,6-Trichlorophenol	5.0	ug/l	1	-	U	Yes
Acenaphthene	1.0	ug/l	1	-	U	Yes
Acenaphthylene	1.0	ug/l	1	-	U	Yes
Acetophenone	2.0	ug/l	1	-	U	Yes
Anthracene	1.0	ug/l	1	-	U	Yes
Atrazine	2.0	ug/l	1	-	U	Yes
Benzaldehyde	5.0	ug/l	1	-	U	Yes
Benzo(a)anthracene	1.0	ug/l	1	-	U	Yes
Benzo(a)pyrene	1.0	ug/l	1	-	U	Yes
Benzo(b)fluoranthene	1.0	ug/l	1	-	U	Yes
Benzo(g,h,i)perylene	1.0	ug/l	1	-	U	Yes
Benzo(k)fluoranthene	1.0	ug/l	1	-	U	Yes
4-Bromophenyl phenyl ether	2.0	ug/l	1	-	U	Yes
Butyl benzyl phthalate	2.0	ug/l	1	-	U	Yes
1,1'-Biphenyl	1.0	ug/l	1	-	U	Yes
2-Chloronaphthalene	2.0	ug/l	1	-	U	Yes
4-Chloroaniline	5.0	ug/l	1	-	U	Yes
Carbazole	1.0	ug/l	1	-	U	Yes
Caprolactam	2.0	ug/l	1	-	U	Yes
Chrysene	1.0	ug/l	1	-	U	Yes
bis(2-Chloroethoxy)methane	2.0	ug/l	1	-	U	Yes
bis(2-Chloroethyl)ether	2.0	ug/l	1	-	U	Yes
bis (2-Chlorois opropyl) ether	2.0	ug/l	1	-	U	Yes

4-Chlorophenyl phenyl ether	2.0	ug/l	1	12	U	Yes
2,4-Dinitrotoluene	1.0	ug/l	1		U	Yes
2,6-Dinitrotoluene	1.0	ug/l	1	2-	U	Yes
3,3'-Dichlorobenzidine	2.0	ug/l	1	-	U	Yes
Dibenzo(a,h)anthracene	1.0	ug/l	1	-	U	Yes
Dibenzofuran	5.0	ug/l	1	-	U	Yes
Di-n-butyl phthalate	2.0	ug/l	1	-	U	Yes
Di-n-octyl phthalate	2.0	ug/l	1		Ų	Yes
Diethyl phthalate	2.0	ug/l	1	_	U	Yes
Dimethyl phthalate	2.0	ug/l	1		U	Yes
bis(2-Ethylhexyl)phthalate	2.0	ug/l	1	-	U	Yes
Fluoranthene	1.0	ug/l	1	-	U	Yes
Fluorene	1.0	ug/l	1	-	U	Yes
Hexachlorobenzene	1.0	ug/l	1	727	U	Yes
Hexachlorobutadiene	1.0	ug/l	1	-	U	Yes
Hexachlorocyclopentadiene	10	ug/l	1	-	UJ	Yes
Hexachloroethane	2.0	ug/l	1	4	U	Yes
Indeno(1,2,3-cd)pyrene	1.0	ug/l	1	-	U	Yes
Isophorone	2.0	ug/l	1	-	U	Yes
1-Methylnaphthalene	1.0	ug/l	1	-	U	Yes
2-Methylnaphthalene	1.0	ug/l	1	(5)	U	Yes
2-Nitroaniline	5.0	ug/l	1	-	U	Yes
3-Nitroaniline	5.0	ug/l	1	-	U	Yes
4-Nitroaniline	5.0	ug/l	1	10.70	U	Yes
Nitrobenzene	2.0	ug/l	1		U	Yes
N-Nitroso-di-n-propylamine	2.0	ug/l	1	•	U	Yes
Nitrosodiphenylamine	5.0	ug/l	1	-	U	Yes
Phenanthrene	1.0	ug/l	1	-	U	Yes
Pyrene	1.0	ug/l	1	(.5)	U	Yes
1,2,4,5-Tetrachlorobenzene	2.0	ug/l	1	-	U	Yes
METHOD:	8270D (SI	M)				
Naphthalene	0.10	ug/l	1	12		Yes
1,4-Dioxane	1.49	ug/l	1	-	-	Yes

Sample location: BMSMC Building 5 Area

Sampling date: 6/17/2016 Matrix: Groundwater

METHOD: 82/0D								
Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable		
2-Chlorophenol	5.1	ug/l	1	-	U	Yes		
4-Chloro-3-methyl phenol	5.1	ug/l	1	-	U	Yes		
2,4-Dichlorophenol	2.0	ug/l	1	-	U	Yes		
2,4-Dimethylphenol	5.1	ug/l	1	•	U	Yes		
2,4-Dinitrophenol	10	ug/l	1	-	U	Yes		
4,6-Dinitro-o-cresol	5.1	ug/l	1	-	U	Yes		
2-Methylphenol	2.0	ug/l	1	-	U	Yes		
3&4-Methylphenol	2.0	ug/l	1	-	U	Yes		
2-Nitrophenol	5.1	ug/l	1	-	Ų	Yes		
4-Nitrophenol	10	ug/l	1	-	U	Yes		
Pentachlorophenol	5.1	ug/l	1	-	U	Yes		
Phenol	2.0	ug/l	1	-	U	Yes		
2,3,4,6-Tetrachlorophenol	5.1	ug/l	1	-	U	Yes		
2,4,5-Trichlorophenol	5.1	ug/l	1	-	Ü	Yes		
2,4,6-Trichlorophenol	5.1	ug/l	1	-	U	Yes		
Acenaphthene	1.0	ug/l	1	-	U	Yes		
Acenaphthylene	1.0	ug/l	1	-	U	Yes		
Acetophenone	2.0	ug/l	1	-	U	Yes		
Anthracene	1.0	ug/l	1	-	U	Yes		
Atrazine	2.0	ug/l	1	-	U	Yes		
Benzaldehyde	5.1	ug/l	1	-	U	Yes		
Benzo(a)anthracene	1.0	ug/l	1	-	U	Yes		
Benzo(a)pyrene	1.0	ug/l	1	-	U	Yes		
Benzo(b)fluoranthene	1.0	ug/l	1	-	U	Yes		
Benzo(g,h,i)perylene	1.0	ug/l	1	-	U	Yes		
Benzo(k)fluoranthene	1.0	ug/l	1	-	U	Yes		
4-Bromophenyl phenyl ether	2.0	ug/l	1	-	U	Yes		
Butyl benzyl phthalate	2.0	ug/i	1	-	U	Yes		
1,1'-Biphenyl	1.0	ug/l	1	-	U	Yes		
2-Chloronaphthalene	2.0	ug/l	1	-	U	Yes		
4-Chloroaniline	1.0	ug/l	1	-	U	Yes		
Carbazole	1.0	ug/l	1	-	U	Yes		
Caprolactam	2.0	ug/l	1	-	U	Yes		
Chrysene	1.0	ug/l	1	-	U	Yes		
bis(2-Chloroethoxy)methane	2.0	ug/l	1	-	U	Yes		
bis(2-Chloroethyl)ether	2.0	ug/l	1	-	U	Yes		
bis (2-Chloro is opropyl) ether	2.0	ug/l	1	-	U	Yes		
4-Chlorophenyl phenyl ether	2.0	ug/l	1	-	U	Yes		

2,4-Dinitrotoluene	1.0	ug/l	1		U	Yes
2,6-Dinitrotoluene	1.0	ug/l	1		U	Yes
3,3'-Dichlorobenzidine	2.0	ug/l	1	-	U	Yes
Dibenzo(a,h)anthracene	1.0	ug/l	1		U	Yes
Dibenzofuran	5.1	ug/l	1		U	Yes
Di-n-butyl phthalate	2.0	ug/l	1	- 2	U	Yes
Di-n-octyl phthalate	2.0	ug/l	1		U	Yes
Diethyl phthalate	2.0	ug/l	1	-	U	Yes
Dimethyl phthalate	2.0	ug/l	1	-	U	Yes
bis(2-Ethylhexyl)phthalate	2.0	ug/l	1	*	U	Yes
Fluoranthene	1.0	ug/l	1	-	U	Yes
Fluorene	1.0	ug/l	1	*)	U	Yes
Hexachlorobenzene	1.0	ug/l	1	-	U	Yes
Hexachlorobutadiene	1.0	ug/l	1	-	U	Yes
Hexachlorocyclopentadiene	10	ug/l	1	5.4	UJ	Yes
Hexachloroethane	2.0	ug/l	1	0.75	U	Yes
Indeno(1,2,3-cd)pyrene	1.0	ug/l	1	_	U	Yes
Isophorone	2.0	ug/l	1	-	U	Yes
1-Methylnaphthalene	1.0	ug/l	1	-	U	Yes
2-Methylnaphthalene	1.0	ug/l	1	-	U	Yes
2-Nitroaniline	5.1	ug/l	1	-	U	Yes
3-Nitroaniline	5.1	ug/l	1	12	U	Yes
4-Nitroaniline	5.1	ug/l	1	-	U	Yes
Nitrobenzene	2.0	ug/l	1	0.75	U	Yes
N-Nitroso-di-n-propylamine	2.0	ug/l	1	-	U	Yes
Nitrosodiphenylamine	5.1	ug/l	1	-	U	Yes
Phenanthrene	1.0	ug/l	1	-	U	Yes
Pyrene	1.0	ug/l	1	-	U	Yes
1,2,4,5-Tetrachlorobenzene	2.0	ug/l	1	6.73	U	Yes
METHOD	9270D (CII	N 4\				
METHOD: Naphthalene	0.10	•	1		U	Yes
1,4-Dioxane	0.334	ug/l	1	7.0	-	Yes
T'4-DIOXALIG	U.334	ug/l	1	_	-	167

Sample location: BMSMC Building 5 Area

Sampling date: 6/17/2016 Matrix: Groundwater

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
2-Chlorophenol	5.0	ug/i	1	-	U	Yes
4-Chloro-3-methyl phenol	5.0	ug/l	1	-	U	Yes
2,4-Dichlorophenol	2.0	ug/l	1	-	U	Yes
2,4-Dimethylphenol	5.0	ug/l	1	-	U	Yes
2,4-Dinitrophenol	10	ug/l	1	-	U	Yes
4,6-Dinitro-o-cresol	5.0	ug/l	1	-	U	Yes
2-Methylphenol	2.0	ug/l	1	-	U	Yes
3&4-Methylphenol	2.0	ug/l	1	-	U	Yes
2-Nitrophenol	5.0	ug/l	1	-	U	Yes
4-Nitrophenol	10	ug/l	1	-	U	Yes
Pentachlorophenol	5.0	ug/l	1	-	U	Yes
Phenol	2.0	ug/l	1	-	U	Yes
2,3,4,6-Tetrachlorophenol	5.0	ug/l	1	-	U	Yes
2,4,5-Trichlorophenol	5.0	ug/l	1	-	U	Yes
2,4,6-Trichlorophenol	5.0	ug/l	1	-	U	Yes
Acenaphthene	1.0	ug/l	1	-	U	Yes
Acenaphthylene	1.0	ug/l	1	-	U	Yes
Acetophenone	2.0	ug/l	1	-	U	Yes
Anthracene	1.0	ug/i	1	•	U	Yes
Atrazine	2.0	ug/l	1	-	U	Yes
Benzaldehyde	5.0	ug/i	1	-	U	Yes
Benzo(a)anthracene	1.0	ug/l	1	-	U	Yes
Benzo(a)pyrene	1.0	ug/l	1	-	U	Yes
Benzo(b)fluoranthene	1.0	ug/l	1	-	U	Yes
Benzo(g,h,i)perylene	1.0	ug/l	1	-	U	Yes
Benzo(k)fluoranthene	1.0	ug/l	1	-	U	Yes
4-Bromophenyl phenyl ether	2.0	ug/l	1	-	U	Yes
Butyl benzyl phthalate	2.0	ug/l	1	-	U	Yes
1,1'-Biphenyl	1.0	ug/l	1	-	U	Yes
2-Chloronaphthalene	2.0	ug/l	1	-	U	Yes
4-Chloroaniline	5.0	ug/l	1	-	U	Yes
Carbazole	1.0	ug/l	1	-	U	Yes
Caprolactam	2.0	ug/l	1	-	U	Yes
Chrysene	1.0	ug/l	1	-	U	Yes
bis(2-Chloroethoxy)methane	2.0	ug/l	1	-	U	Yes
bis(2-Chloroethyl)ether	2.0	ug/l	1	-	U	Yes
bis (2-Chloro is opropyl) ether	2.0	ug/l	1	-	U	Yes
4-Chlorophenyl phenyl ether	2.0	ug/l	1	-	U	Yes

2,4-Dinitrotoluene	1.0	ug/l	1	-	U	Yes
2,6-Dinitrotoluene	1.0	ug/l	1		U	Yes
3,3'-Dichlorobenzidine	2.0	ug/l	1	-	U	Yes
Dibenzo(a,h)anthracene	1.0	ug/l	1		U	Yes
Dibenzofuran	5.0	ug/l	1	-	U	Yes
Di-n-butyl phthalate	2.0	ug/l	1	-	U	Yes
Di-n-octyl phthalate	2.0	ug/l	1	-	U	Yes
Diethyl phthalate	2.0	ug/l	1		U	Yes
Dimethyl phthalate	2.0	ug/l	1	-	U	Yes
bis(2-Ethylhexyl)phthalate	2.0	ug/l	1	-	U	Yes
Fluoranthene	1.0	ug/l	1	-	U	Yes
Fluorene	1.0	ug/l	1	-	Ų	Yes
Hexachlorobenzene	1.0	ug/l	1		U	Yes
Hexachlorobutadiene	1.0	ug/l	1		U	Yes
Hexachlorocyclopentadiene	10	ug/l	1	-	UJ	Yes
Hexachloroethane	2.0	ug/l	1		U	Yes
Indeno(1,2,3-cd)pyrene	1.0	ug/l	1	-	U	Yes
Isophorone	2.0	ug/l	1	4. 5 0	U	Yes
1-Methylnaphthalene	1.0	ug/l	1	-	U	Yes
2-Methylnaphthalene	1.0	ug/l	1	-	U	Yes
2-Nitroaniline	5.0	ug/l	1	-	U	Yes
3-Nitroaniline	5.0	ug/l	1	-	U	Yes
4-Nitroaniline	5.0	ug/l	1	-	U	Yes
Nitrobenzene	2.0	ug/l	1		U	Yes
N-Nitroso-di-n-propylamine	2.0	ug/l	1	-	U	Yes
Nitrosodiphenylamine	5.0	ug/l	1	-	U	Yes
Phenanthrene	1.0	ug/l	1	-	U	Yes
Pyrene	1.0	ug/l	1	-	U	Yes
1,2,4,5-Tetrachlorobenzene	2.0	ug/l	1		U	Yes
METHOD:	8270D (SII	\A\				
Naphthalene	0.10	ug/l	1	-	U	Yes
1,4-Dioxane	0.398	ug/i	i	-	-	Yes
I, T DIOXOIIC	0.556	ug/1	*			163

Sample location: BMSMC Building 5 Area

Sampling date: 6/20/2016 Matrix: Groundwater

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
2-Chlorophenol	6.0	ug/l	1	-	U	Yes
4-Chloro-3-methyl phenol	6.0	ug/l	1	-	U	Yes
2,4-Dichlorophenol	2.4	ug/l	1	-	U	Yes
2,4-Dimethylphenol	6.0	ug/l	1	-	U	Yes
2,4-Dinitrophenol	12	ug/l	1	-	U	Yes
4,6-Dinitro-o-cresol	6.0	ug/l	1	-	U	Yes
2-Methylphenol	2.4	ug/l	1	-	U	Yes
3&4-Methylphenol	2.4	ug/l	1	-	U	Yes
2-Nitrophenol	6.0	ug/l	1	-	U	Yes
4-Nitrophenol	12	ug/l	1	-	U	Yes
Pentachlorophenol	6.0	ug/l	1	-	U	Yes
Phenol	2.4	ug/l	1	-	U	Yes
2,3,4,6-Tetrachlorophenol	6.0	ug/i	1	-	U	Yes
2,4,5-Trichlorophenol	6.0	ug/l	1	-	U	Yes
2,4,6-Trichlorophenol	6.0	ug/l	1	-	U	Yes
Acenaphthene	1.2	ug/l	1	-	U	Yes
Acenaphthylene	1.2	ug/l	1	-	U	Yes
Acetophenone	2.4	ug/l	1	-	U	Yes
Anthracene	1.2	ug/l	1	-	U	Yes
Atrazine	2.4	ug/l	1	-	U	Yes
Benzaldehyde	6.0	ug/l	1	-	U	Yes
Benzo(a)anthracene	1.2	ug/l	1	-	U	Yes
Benzo(a)pyrene	1.2	ug/l	1	-	U	Yes
Benzo(b)fluoranthene	1.2	ug/l	1	-	U	Yes
Benzo(g,h,i)perylene	1.2	ug/l	1	-	U	Yes
Benzo(k)fluoranthene	1.2	ug/l	1	-	U	Yes
4-Bromophenyl phenyl ether	2.4	ug/l	1	-	U	Yes
Butyl benzyl phthalate	2.4	ug/l	1	-	U	Yes
1,1'-Biphenyl	1.2	ug/l	1	-	U	Yes
2-Chloronaphthalene	2.4	ug/l	1	-	U	Yes
4-Chloroaniline	6.0	ug/l	1	-	U	Yes
Carbazole	1.2	ug/l	1	-	U	Yes
Caprolactam	2.4	ug/l	1	-	U	Yes
Chrysene	1.2	ug/l	1	-	U	Yes
bis(2-Chloroethoxy)methane	2.4	ug/l	1	-	U	Yes
bis(2-Chloroethyl)ether	2.4	ug/l	1	-	U	Yes
bis (2-Chloro is opropyl) ether	2.4	ug/l	1	-	U	Yes
4-Chlorophenyl phenyl ether	2.4	ug/l	1	-	U	Yes

2,4-Dinitrotoluene	1.2	ug/l	1	-	υ	Yes
2,6-Dinitrotoluene	1.2	ug/l	1	-	U	Yes
3,3'-Dichlorobenzidine	2.4	ug/l	1	-	U	Yes
Dibenzo(a,h)anthracene	1.2	ug/l	1	147	U	Yes
Dibenzofuran	6.0	ug/l	1	-	U	Yes
Di-n-butyl phthalate	2.4	ug/l	1	-	U	Yes
Di-n-octyl phthalate	2.4	ug/l	1	-	U	Yes
Diethyl phthalate	2.4	ug/l	1	57.	U	Yes
Dimethyl phthalate	2.4	ug/l	1	-	U	Yes
bis (2-Ethylhexyl) phthalate	2.4	ug/l	1	-	U	Yes
Fluoranthene	1.2	ug/l	1	-	U	Yes
Fluorene	1.2	ug/l	1		U	Yes
Hexachlorobenzene	1.2	ug/l	1	-	U	Yes
Hexachlorobutadiene	1.2	ug/l	1	_	U	Yes
Hexachlorocyclopentadiene	12	ug/l	1	-	UJ	Yes
Hexachloroethane	2.4	ug/l	1	57.5	U	Yes
Indeno(1,2,3-cd)pyrene	1.2	ug/l	1	-	U	Yes
Isophorone	2.4	ug/l	1	-	U	Yes
1-Methylnaphthalene	1.2	ug/l	1	-	U	Yes
2-Methylnaphthalene	1.2	ug/l	1	-	U	Yes
2-Nitroaniline	6.0	ug/l	1		U	Yes
3-Nitroaniline	6.0	ug/l	1	227	U	Yes
4-Nitroaniline	6.0	ug/l	1	-	U	Yes
Nitrobenzene	2.4	ug/l	1	7	U	Yes
N-Nitroso-di-n-propylamine	2.4	ug/l	1	-	U	Yes
Nitrosodiphenylamine	6.0	ug/l	1	-	U	Yes
Phenanthrene	1.2	ug/l	1	-	U	Yes
Pyrene	1.2	ug/l	1	-	U	Yes
1,2,4,5-Tetrachlorobenzene	2.4	ug/l	1	-	U	Yes
METHOD:	8270D (SI	M				
Naphthalene	0.12	ug/l	1	-	U	Yes
1,4-Dioxane	0.12	ug/l	10	-	Ü	Yes
wy - writerial its		-6/ ·	20		~	,

Sample location: BMSMC Building 5 Area

Sampling date: 6/17/2016 Matrix: Groundwater

WETHOD:			_			
Analyte Name	Result		Dilution Factor	Lab Flag		=
2-Chlorophenol	5.5	ug/l	1	-	U	Yes
4-Chloro-3-methyl phenol	5.5	ug/l	1	-	U	Yes
2,4-Dichlorophenol	2.2	ug/l	1	-	U	Yes
2,4-Dimethylphenol	5.5	ug/l	1	-	U	Yes
2,4-Dinitrophenol	11	ug/l	1	-	Ų	Yes
4,6-Dinitro-o-cresol	5.5	ug/l	1	-	U	Yes
2-Methylphenol	2.2	ug/l	1	-	U	Yes
3&4-Methylphenol	2.2	ug/l	1	-	U	Yes
2-Nitrophenol	5.5	ug/l	1	-	U	Yes
4-Nitrophenol	11	ug/l	1	-	U	Yes
Pentachlorophenol	5.5	ug/l	1	-	U	Yes
Phenol	2.2	ug/l	1	-	U	Yes
2,3,4,6-Tetrachlorophenol	5.5	ug/l	1	-	U	Yes
2,4,5-Trichlorophenol	5.5	ug/l	1	-	U	Yes
2,4,6-Trichlorophenol	5.5	ug/l	1	-	U	Yes
Acenaphthene	1.1	ug/l	1	-	U	Yes
Acenaphthylene	1.1	ug/l	1	-	U	Yes
Acetophenone	2.2	ug/l	1	-	U	Yes
Anthracene	1.1	ug/l	1	-	U	Yes
Atrazine	2.2	ug/l	1	-	U	Yes
Benzaldehyde	5.5	ug/l	1	-	U	Yes
Benzo(a)anthracene	1.1	ug/l	1	-	U	Yes
Benzo(a)pyrene	1.1	ug/l	1	-	U	Yes
Benzo(b)fluoranthene	1.1	ug/l	1	-	U	Yes
Benzo(g,h,i)perylene	1.1	ug/l	1	-	U	Yes
Benzo(k)fluoranthene	1.1	ug/l	1	-	U	Yes
4-Bromophenyl phenyl ether	2.2	ug/l	1	-	U	Yes
Butyl benzyl phthalate	2.2	ug/l	1	-	U	Yes
1,1'-Biphenyl	1.1	ug/l	1	-	U	Yes
2-Chloronaphthalene	2.2	ug/l	1	-	U	Yes
4-Chloroaniline	5.5	ug/l	1	-	U	Yes
Carbazole	1.1	ug/l	1	-	U	Yes
Caprolactam	2.2	ug/l	1	-	U	Yes
Chrysene	1.1	ug/l	1	-	U	Yes
bis(2-Chloroethoxy)methane	2.2	ug/l	1	-	U	Yes
bis(2-Chloroethyl)ether	2.2	ug/l	1	-	U	Yes
bis(2-Chloroisopropyl)ether	2.2	ug/l	1	_	U	Yes
4-Chlorophenyl phenyl ether	2.2	ug/l	1	-	U	Yes
		_				

2,4-Dinitrotoluene	1.1	ug/l	1	-	U	Yes
2,6-Dinitrotoluene	1.1	ug/l	1	7	U	Yes
3,3'-Dichlorobenzidine	2.0	ug/l	1	-	U	Yes
Dibenzo(a,h)anthracene	1.1	ug/l	1	-	U	Yes
Dibenzofuran	5.5	ug/l	1		U	Yes
Di-n-butyl phthalate	2.2	ug/l	1	-	U	Yes
Di-n-octyl phthalate	2.2	ug/l	1	-	U	Yes
Diethyl phthalate	2.2	ug/l	1	-	Ų	Yes
Dimethyl phthalate	2.2	ug/l	1	-	U	Yes
bis(2-Ethylhexyl)phthalate	2.2	ug/l	1		U	Yes
Fluoranthene	1.1	ug/l	1		U	Yes
Fluorene	1.1	ug/l	1		U	Yes
Hexachlorobenzene	1.1	ug/l	1	7.5	U	Yes
Hexachlorobutadiene	1.1	ug/l	1	-	U	Yes
Hexachlorocyclopentadiene	10	ug/l	1		IJ	Yes
Hexachloroethane	2.2	ug/l	1	-	U	Yes
Indeno(1,2,3-cd)pyrene	1.1	ug/l	1	-	U	Yes
Isophorone	2.0	ug/i	1	17.0	U	Yes
1-Methylnaphthalene	1.1	ug/l	1		U	Yes
2-Methylnaphthalene	1.1	ug/l	1	-	U	Yes
2-Nitroaniline	5.5	ug/l	1	35.	U	Yes
3-Nitroaniline	5.5	ug/l	1	-	Ų	Yes
4-Nitroaniline	5.5	ug/l	1	.*	Ų	Yes
Nitrobenzene	2.2	ug/l	1	-	U	Yes
N-Nitroso-di-n-propylamine	2.2	ug/l	1	-	U	Yes
Nitrosodiphenylamine	5.5	ug/l	1	77	U	Yes
Phenanthrene	1.1	ug/l	1	-	U	Yes
Pyrene	1.1	ug/l	1		U	Yes
1,2,4,5-Tetrachlorobenzene	2.2	ug/l	1	-	U	Yes

METHOD: 8270D (SIM)

Naphthalene	0.11	ug/l	1	-	U	Yes
1,4-Dioxane	0.646	ug/l	1	-	-	Yes

	Project Number:_JC22613
	Date: June_16-20,_2016
	Shipping Date:June_16,_2016
	EPA Region:2_
REVIEW OF SEMIVOLATION	LE ORGANIC PACKAGE
The following guidelines for evaluating volatile validation actions. This document will assist the make more informed decision and in better ser results were assessed according to USEPA following order of precedence: EPA Hazardou 2015 –Revision 0. Semivolatile Data Validation. The on the data review worksheets are from the proted.	ne reviewer in using professional judgment to ving the needs of the data users. The sample data validation guidance documents in the is Waste Support Section, SOP HW-35A, July e QC criteria and data validation actions listed
The hardcopied (laboratory name) _Accutest reviewed and the quality control and performance included:	data package received has been data summarized. The data review for SVOCs
Lab. Project/SDG No.:JC22613 No. of Samples:8_Full_scan/8_SIM	
Trip blank No.:	
Field blank No.:	
Equipment blank No.:	
Field duplicate No.:	
X Data Completeness	X Laboratory Control Spikes
X Holding Times	X Field Duplicates
X GC/MS Tuning	X Calibrations
X Internal Standard Performance	X Compound Identifications
X Blanks	X Compound Quantitation
X Surrogate RecoveriesX Matrix Spike/Matrix Spike Duplicate	X Quantitation Limits
Overall Comments:_ABN_TCL_list_by_method_SW8 _analyzed_by_method_SW846-8270D_(SIM)	
Definition of Qualifiers:	
J- Estimated results	
U- Compound not detected	
R- Rejected data // /	
UJ- Estimated nongetect/	
Reviewer: (afail Mant	
Date:July_21,_2016	

DATA COMPLETENESS

MISSING INFORMATION	DATE LAB. CONTACTED	DATE RECEIVED
- 2000	1000	
·····		
		**** **** **** *****

All criteria were met _	X_	_
Criteria were not met		
and/or see below		

HOLDING TIMES

The objective of this parameter is to ascertain the validity of the results based on the holding time of the sample from time of collection to the time of analysis.

Complete table for all samples and note the analysis and/or preservation not within criteria

SAMPLE ID	DATE SAMPLED	DATE EXTRACTED/ANALYZED	рН	ACTION
All samples extr	acted and ana	alyzed within method recom	men	ded holding time. Samples properly

Cooler temperature ((Criteria: 4 <u>+</u> 2 °C)):3.9°C
----------------------	-----------------------------	---------

Actions

Results will be qualified based on the criteria of the following Table:

Table 1. Holding Time Actions for Semivolatile Analyses

Table 1. Holding Time Actions for Semivolatile Analyses Action				
Matrix	Preserved	Criteria	Detected Associated Compounds	Non-Detected Associated Compounds
	No	≤ 7 days (for extraction) ≤ 40 days (for analysis)	Use professi	onal judgment
	No	> 7 days (for extraction) > 40 days (for analysis)	J	Use professional judgment
Aqueous	Yes	≤ 7 days (for extraction) ≤ 40 days (for analysis)	No qualification	
	Yes	> 7 days (for extraction) > 40 days (for analysis)	J	UJ
	Yes/No	Grossly Exceeded	J	UJ or R
	No	≤ 14 days (for extraction) ≤ 40 days (for analysis)	Use profession	onal judgment
Non-Aqueous	No	> 14 days (for extraction) > 40 days (for analysis)	I.	Use professional judgment
	Yes	≤ 14 days (for extraction) ≤ 40 days (for analysis)	No qualification	
	Yes	> 14 days (for extraction) > 40 days (for analysis)	J	ÜJ
	Yes/No	Grossly Exceeded	J	ÚJ or R

All criteria were met _X
Criteria were not met see below

GC/MS TUNING

The assessment of the tuning results is to determine if the sample instrumentation is within the standard tuning QC limits

_X__ The DFTPP performance results were reviewed and found to be within the specified criteria.

_X__ DFTPP tuning was performed for every 12 hours of sample analysis.

If no, use professional judgment to determine whether the associated data should be accepted, qualified or rejected.

Notes: These requirements do not apply when samples are analyzed by the Selected Ion Monitoring (SIM) technique.

All mass spectrometer conditions must be identical to those used during the sample analysis. Background subtraction actions resulting in spectral distortion are unacceptable

Notes: No data should be qualified based of DFTPP failure.

The requirement to analyze the instrument performance check solution is optional when analysis of PAHs/pentachlorophenol is to be performed by the SIM technique.

LIST	tne	samples	апестеа:
	<u> </u>		

Actions:

- If sample are analyzed without a preceding valid instrument performance check or are analyzed 12 hours after the Instrument Performance Check, qualify all data in those samples as unusable (R).
- 2. If ion abundance criteria are not met, use professional judgment to determine to what extent the data may be utilized.
- 3. State in the Data Review Narrative, decisions to use analytical data associated with DFTPP instrument performance checks not meeting the contract requirements.
- 4. Use professional judgment to determine if associated data should be qualified based on the spectrum of the mass calibration compounds.

All criteria were met _	_X
Criteria were not met	
and/or see below	

INITIAL CALIBRATION VERIFICATION

Compliance requirements for satisfactory instrument calibration are established to ensure that the instrument is capable of producing and maintaining acceptable quantitative data.

	06/22/16_(Scan) GCMS3E Aqueous/low
Date of initial calibration: Instrument ID numbers: Matrix/Level:	GCMSM
	06/20/16_(SIM) GCMS3M _Aqueous/low

	DATE	LAB ID#	FILE	CRITERIA OUT RFs, %RSD, %D, r	COMPOUND	SAMPLES AFFECTED			
200	Initial and initial calibration verification meets the method and guidance validation document performance criteria. Other instrument used for the analysis of QC samples for this job. The QC samples analyzed were not validated								

Actions:

Qualify the initial calibration analytes listed in Table 2 using the following criteria:

Table 3. Initial Calibration Actions for Semivolatile Analysis

		Action
Criteria	Detect	Non-detect
Initial Calibration not performed at specified frequency and sequence	Use professional judgment R	Use professional judgment R
Initial Calibration not performed at the specified concentrations	J	ÜJ
RRF < Minimum RRF in Table 2 for target analyte	Use professional judgment J J or R	R
RRF > Minimum RRF in Table 2 for target analyte	No qualification	No qualification
%RSD > Maximum %RSD in Table 2 for target analyte	Ţ	Use professional judgment
%RSD ≤ Maximum %RSD in Table 2 for target analyte	No qualification	No qualification

Initial Calibration

Table 2. RRF, %RSD, and %D Acceptance Criteria in Initial Calibration and CCV for Semivolatile Analysis

Analyte	Minimum RRF	Maximum %RSD	Opening Maximum %D ¹	Opening Maximum %D ¹
1,4-Dioxane	0.010	40,0	-40.0	- 50.0
Benzaldehyde	0.100	40,0	= 40.0	= 50.0
Phenol	0.080	20.0	= 20.0	= 25.0
Bis(2-chloroethyl)ether	0.100	20.0	= 20.0	= 25.0
2-Chlorophenol	0,200	20,0	=20,0	= 25.0
2-Methylphenol	0,010	20.0	= 20.0	= 25.0
3-Methylphenol	0.010	20.0	-20.0	±25.0
2,2'-Oxybis-(1-chloropropane)	0.010	20.0	= 25.0	= 50.0
Acetophenone	0.060	20,0	= 20.0	= 25,0
4-Methylphenol	0.010	20.0	- 20.0	-25.0
N-Nitroso-di-n-propylamine	0.080	20.0	- 25.0	=25.0
Hexachloroethane	0.100	20.0	-20.0	-25.0
Nitrobenzene	0.090	20.0	= 20.0	±25.0
Isophorone	0,100	20.0	= 20.0	- 25.0
2-Nitrophenol	0.060	20.0	- 20.0	- 25.0
2,4-Dimethylphenol	0.050	20.0	25.0	= 50.0
Bis(2-chloroethoxy)methane	0.080	20.0	-20.0	-25.0
2,4-Dichlorophenol	0.060	20.0	= 20.0	= 25.0
Naphthalene	0.200	20.0	20.0	- 25.0
4-Chloroaniline	0.010	40.0	±40.0	± 50.0
Hexachlorobutadiene	0.040	20,0	= 20.0	= 25.0
Caprolactam	0.010	40.0	= 30.0	- 50.0
4-Chloro-3-methy3phenol	0.040	20.0	= 20.0	= 25.0
2-MethyInaphthalene	0.100	20,0	-20.0	±25.0
Hexachlorocyclopentadiene	0.010	40.0	±40.0	= 50.0
2,4,6-Trichtorophenol	0.090	20.0	-20.0	- 25.0
2,4,5-Trichlorophenol	0,100	20,0	= 20.0	= 25.0
1,1'-Biphenyl	0.200	20,0	= 20.0	- 25.0

Analyte	Minimum RRF	Maximum %RSD	Opening Maximum %D ¹	Opening Maximum %D ¹
2-C'hloronaphthalene	0.300	20.0	-20.0	25.0
2-Nitroaniline	0.060	20.0	-25.0	± 25.0
Dimethylphthalate	0.300	20.0	= 25.0	±25.0
2,6-Dinitrotoluene	0.080	20.0	= 20.0	-25.0
Acenaphthylene	0.400	20.0	= 20.0	= 25.0
3-Nitroaniline	0.010	20.0	-25.0	-50.0
Acenaphthene	0.200	20.0	- 20.0	- 25.0
2.4-Dinitrophenol	0.010	40.0	- 50.0	-50.0
4-Nitrophenol	0.010	40.0	= 40.0	= 50.0
Dibenzofuran	0.300	20.0	= 20.0	-25.0
2,4-Dinitrotoluene	0.070	20,0	-20.0	- 25.0
Diethylphthalate	0.300	20.0	= 20.0	- 25.0
1,2,4,5-Tetrachlorobenzene	0.100	20.0	= 20.0	-25.0
4-Chlorophenyl-phenylether	0.100	20,0	- 20.0	- 25.0
Fluorene	0.200	20.0	= 20.0	- 25.0
4-Nitroaniline	0,010	40.0	-40.0	-50.0
4,6-Dinitro-2-methylphenol	0.010	40.0	=30.0	-50.0
4-Bromophenyl-phenyl ether	0.070	20.0	= 20.0	± 25.0
N-Nitrosodiphenylamine	0.100	20.0	- 20.0	- 25.0
Hexachlorobenzene	0.050	20,0	-20.0	- 25.0
Atrazine	0.010	40.0	-25.0	-50.0
Pentachlorophenol	0.010	40.0	-40.0	-50.0
Phenanthrene	0.200	20.0	-20.0	-25.0
Anthracene	0,200	20.0	- 20.0	- 25.0
Carbazole	0.050	20.0	- 20.0	-25.0
Di-n-butylphthalate	0.500	20.0	- 20,0	-25.0
Fluoranthene	0.100	20,0	- 20.0	-25.0
Pyrene	0,400	20.0	-25.0	-50.0
Butylbenzylphthalate	0,100	20.0	-25.0	- 50.0

Analyte	Minimum RRF	Maximum %RSD	Opening Maximum %D ¹	Opening Maximum %D ¹
3,3'-Dichlorobenzidine	0.010	40.0	- 40.0	- 50.0
Benzo(a)anthracene	0.300	20.0	- 20.0	- 25.0
Chrysene	0.200	20.0	= 20.0	- 50.0
Bis(2-ethylhexyl) phthalate	0,200	20.0	-25.0	-50.0
Di-n-octylphthalate	0.010	40,0	=40.0	- 50.0
Benzo(b)fluoranthene	0.010	20.0	± 25.0	= 50.0
Benzo(k)fluoranthene	0.010	20,0	=25.0	- 50.0
Benzo(a)pyrene	0.010	20.0	- 20.0	- 50.0
Indeno(1,2,3-cd)pyrene	0.010	20.0	= 25.0	 50.0
Dibenzo(a,h)anthracene	0.010	20.0	± 25.0	- 50.0
Benzo(g,h,i)perylene	0.010	20.0	-30.0	- 50.0
2,3,4,6-Tetrachlorophenol	0.040	20.0	= 20.0	= 50.0
Naphthalene	0,600	20.0	÷25.0	- 25.0
2-Methylnaphthalene	0.300	20.0	± 20.0	- 25.0
Acenaphthylene	0.900	20.0	± 20.0	- 25.0
Acenaphthene	0,500	20,0	- 20.0	- 25.0
Fluorene	0.700	20,0	= 25.0	= 50.0
Phenanthrene	0,300	20.0	- 25.0	= 50.0
Anthracene	0.400	20.0	± 25.0	= 50,0
Fluoranthene	0.400	20.0	- 25.0	= 50.0
Pyrene	0.500	20.0	=30.0	= 50.0
Benzo(a)anthracene	0,400	20.0	=25.0	= 50.0
Chyrsene	0,400	20.0	= 25.0	= 50.0
Benzo(b)fluoranthene	0,100	20.0	=,30.0	- 50.0
Benzo(k)fluoranthene	0.100	20.0	-30.0	= 50.0
Benzo(a)pyrene	0.100	20.0	= 25.0	= 50.0
Indeno(1,2,3-cd)pyrene	0.100	20.0	-40.0	= 50.0
Dibenzo(a,h)anthracene	010.0	25.0	-40.0	- 50.0
Benzo(g,h,i)perylene	0.020	25.0	-40.0	- 50.0

Pentachlorophenol	010.0	40,0	- 50.0	- 50.0	٦		
Deuterated Monitoring Compounds							

Analyte	Minimum RRF	Maximum %RSD	Opening Maximum %D ¹	Closing Maximum %D
1,4-Dioxane-d _x	0.010	20.0	= 25.0	-50.0
Phenol-d ₃	0.010	20.0	- 25.0	+25.0
Bis-(2-chloroethyl)ether-d _k	0,100	20.0	- 20_0	-25.0
2-Chlorophenol-di	0.200	20.0	-20.0	-25.0
4-Methylphenol-d _x	0.010	20.0	-20.0	- 25.0
4-Chloroaniline-d ₄	0.010	40.0	-40.0	- 50.0
Nitrobenzene-d ₃	0.050	20.0	~20.0	± 25.0
2-Nitrophenol-d ₁	0.050	20.0	-20.0	- 25.0
2,4-Dichlorophenol-d:	0.060	20.0	-20,0	+25.0
Dimethylphthalate-d ₆	0,300	20.0	-20.0	-25.0
Acenaphthylene-d ₈	0,400	20.0	- 20.0	-25.0
4-Nitrophenol-d ₁	0.010	40.0	-40.0	- 50.0
Fluorene-d ₁₀	0,100	20.0	= 20.0	=25.0
4,6-Dinitro-2-methylphenol-d2	0.010	40.0	-30,0	- 50.0
Anthracene-d ₁₀	0,300	20.0	± 20.0	-25.0
Pyrene-d ₁₀	0.300	20.0	-25.0	- 50.0
Benzo(a)pyrene-d ₁₂	0.010	20.0	- 20.0	- 50.0
Fluoranthene-d ₁₀ (SIM)	0.400	20.0	± 25.0	-50.0
2-Methylnaphthalene-d ₁₀ (SIM)	0.300	20.0	- 20.0	-25.0

If a closing CCV is acting as an opening CCV, all target analytes must meet the requirements for an opening CCV.

Note: If analysis by SIM technique is requested for PAH/pentachlorophenols, calibration standards analyzed at 0.10, 0.20, 0.40, 0.80, and 1.0 ng/uL for each target compound of interest and the associated DMCs. Pentachlorophenol will require only a four point initial calibration at 0.20, 0.40, 0.80, and 1.0 ng/uL.

All criteria were met	_
Criteria were not met	
and/or see belowX	_

CONTINUING CALIBRATION VERIFICATION

Compliance requirements for satisfactory instrument calibration are established to ensure that the instrument is capable of producing and maintaining acceptable quantitative data.

Date of initial calibration:	06/22/16_(Scan)
Date of initial calibration verification	ion (ICV):_06/22/16
Date of continuing calibration ver	ification (CCV):_07/01/16
Date of closing CCV:	
Instrument ID numbers:	GCMS3E
	Aqueous/low
Date of initial calibration:	06/20/16_(SIM)
Date of initial calibration verificati	on (ICV):_06/20-21/16
	ification (CCV):_06/26/16;_07/01/16_
Instrument ID numbers:	GCMS3M
	Aqueous/low
Date of initial calibration:	06/16/16_(Scan)
	on (ICV):_06/16-17/16;_06/19/16
	ification (CCV):_06/23/16
	-
Instrument ID numbers:	GCMSM
	Aqueous/low

DATE	LAB FILE ID#	CRITERIA OUT RFs, %RSD, %D, r	COMPOUND	SAMPLES AFFECTED
GCMS3E				
07/01/16	cc3699-25	-41.5	Hexachlorocyclopentadiene	JC22613-2 to -&
GCMSM				
6/23/16	cc5307-50	24.9	1,4-dioxane*	JC22613-1
		-24.2	2,3,4,6-tetrachlorophenoi	
	cc5308-50	-24.1	Atrazine*	

Note: Initial and continuing calibration verifications meet the method and guidance document required performance criteria except in the cases described in this document. Analytes not meeting the continuing calibration verification method performance criteria and validation guidance document performance criteria qualified as estimated (J) or (UJ) in affected samples.

No closing calibration verification included in data package. No action taken, professional judgment.

^{*} Analytes not meeting the continuing calibration verification method performance criteria but were within the validation guidance document performance criteria. No action taken.

Actions:

Notes: Verify that the CCV is run at the required frequency (an opening and closing CCV must be run within 12-hour period).

All DMCs must meet the RRF values given in Table 2. No qualification of the data is necessary on DMCs RRF and %RSD/%D alone. Use professional judgment to evaluate DMCs and %RSD/%D data in conjunction with DMCs recoveries to determine the need for qualification of the data.

Qualify the initial calibration analytes listed in Table 2 using the following criteria in the CCVs:

Table 4. CCV Actions for Semivolatile Analysis

Criteria for Opening CCV	Criteria for Closing CCV -	Ac	tion
Criteria for Opening C.C.V	Criteria for Closing CCV	Detect	Non-detect
CCV not performed at required frequency and sequence	CCV not performed at required frequency	Use professional judgment R	Use professional judgment R
CCV not performed at specified concentration	CCV not performed at specified concentration	Use professional judgment	Use professional judgment
RRF < Minimum RRF in Table 2 for target analyte	RRF < Minimum RRF in Table 2 for target analyte	Use professional judgment J or R	R
RRF ≥ Minimum RRF in Table 2 for target analyte	RRF > Minimum RRF in Table 2 for target analyte	No qualification	No qualification
%D outside the Opening Maximum %D limits in Table 2 for target analyte	%D outside the Closing Maximum %D limits in Table 2 for target analyte	ĵ	UJ
%D within the inclusive Opening Maximum %D limits in Table 2 for target analyte	%D within the inclusive Closing Maximum %D limits in Table 2 for target analyte	No qualification	No qualification

All criteria were met _X_	
Criteria were not met	
and/or see below	_

BLANK ANALYSIS RESULTS (Sections 1 & 2)

The assessment of the blank analysis results is to determine the existence and magnitude of contamination problems. The criteria for evaluation of blanks apply only to blanks associated with the samples, including trip, equipment, and laboratory blanks. If problems with any blanks exist, all data associated with the case must be carefully evaluated to determine whether or not there is an inherent variability in the data for the case, or if the problem is an isolated occurrence not affecting other data.

List the contamination in the blanks below. High and low levels blanks must be treated separately.

Notes: The concentration of non-target compounds in all blanks must be less than or equal to 10 ug/L.

The concentration of target compounds in all blanks must be less than its CRQL listed in the method.

Samples taken from a drinking water tap do not have and associated field blank.

Laboratory blanks

DATE Analyzed	LAB ID	LEVEL! MATRIX	COMPOUND	CONCENTRATION UNITS
_No_target_ana	lytes_detected	_in_method_bla	inks	
.				
Field/Equipment	t/Trip blank			
DATE ANALYZED	LAB ID	LEVEL/ MATRIX	COMPOUND	CONCENTRATION UNITS
_No_field/trip/ed	quipment_blank	s_analyzed_wit	h_this_data_package	
2.7765				
			300	
				01 - 3 - 100 - 100 - 100 - 100 - 100 - 100 - 100 - 100 - 100 - 100 - 100 - 100 - 100 - 100 - 100 - 100 - 100 -

All criteria were met _X_	
Criteria were not met	
and/or see below	

BLANK ANALYSIS RESULTS (Section 3)

Blank Actions

Qualify samples based on the criteria summarized in Table 5:

Table 5. Blank and TCLP/SPEP LEB Actions for Semivolatile Analysis

Blank Type	Blank Result	Sample Result	Action
	Detect	Non-detect	No qualification
	< CRQL	< CRQL	Report at CRQL and qualify as non-detect (U)
		> CRQL	Use professional judgment
		< CRQL	Report at CRQL and qualify as non-detect (U)
Method,	•	> CRQL but < Blank Result	Report at sample results and qualify as non-detect (U) or as unusable (R)
TCLP/SPLP LEB, Field		≥ CRQL and ≥ Blank Result	Use professional judgment
	Grossly high	Detect	Report at sample results and qualify as unusable (R)
	TIC > 5.0 ug/L (water) or 0.0050 mg/L (TCLP leachate) or TIC > 170 ug/Kg (soil)	Detect	Use professional judgment

List samples qualified

CONTAMINATION SOURCE/LEVEL	COMPOUND	CONC/UNITS	AL/UNITS	SQL	AFFECTED SAMPLES

All criteria were met	X
Criteria were not met	ST. W
and/or see below	

SURROGATE SPIKE RECOVERIES - DEUTERATED MONITORING COMPOUNDS (DMCs)

Laboratory performance of individual samples is established by evaluation of surrogate spike recoveries – deuterated monitoring compounds. All samples are spiked with surrogate compounds prior to sample analysis. The accuracy of the analysis is measured by the surrogate percent recovery. Since the effects of the sample matrix are frequently outside the control of the laboratory and may present relatively unique problems, the validation of data is frequently subjective and demands analytical experience and professional judgment.

Notes: Recoveries for DMCs in samples and blanks must be within the limits specified in Table 6.

The recovery limits for any of the compounds listed in Table 6 may be expanded at any time during the period of performance if USEPA determines that the limits are too restrictive.

If a DMC is not added in the samples and blanks or the concentrations of DMCs in the samples and blank not the specified, use professional judgment in qualifying the data.

Table 7. DMC Actions for Semivolatile Analysis

Circuit	Action		
Criteria	Detect	Non-detect	
%R < 10% (excluding DMCs with 10% as a lower acceptance limit)	J-	R	
10% ≤ %R (excluding DMC's with 10% as a lower acceptance limit) ≤ Lower Acceptance Limit	J÷	ເສ	
Lower Acceptance limit < %R < Upper Acceptance Limit	No qualification	No qualification	
%R > Upper Acceptance Limit	J)	No qualification	

List the percent recoveries (%Rs) which do not meet the criteria for DMCs (surrogate) recovery.

Matrix:___Groundwater_____

SAMPLE ID SURROGATE COMPOUND ACTION

_DMCs_meet_the_required_criteria._Non-deuterated_surrogates_added_to_the_samples_were__
_within_laboratory_recovery_limits.______

Table 8. Semivolatile DMCs and the Associated Target Analytes

1,4-Dioxane-d ₈ (DMC-1)	Phenot-ds (DMC-2)	Bis(2-Chloroethyl) ether-d ₈
1,4-Dioxane	December 2	(DMC-3)
A.4-170Nane	Benzaldehyde Phenol	Bis(2-chloroethyl)ether
	rnenoi	2,2'-Oxybis(1-chloropropane)
		Bis(2-chloroethoxy)methane
2-Chlorophenol-d4(DMC-4)	4-Methylphenol-d ₈ (DMC-5)	4-Chloroaniline-d4(DMC-6)
2-Chlorophenol	2-Methylphenol	4-Chloroaniline
	3-Methylphenol	Hexachlorocyclopentadiene
	4-Methylphenol	Dichlorobenzidine
	2,4-Dimethylphenol	i
Nitrobenzene-d ₅ (DMC-7)	2-Nitrophenol-d ₄ (DMC-8)	2,4-Dichforophenol-d3(DMC-9)
Acetophenone	Isophorone	2,4-Dichlorophenol
N-Nitroso-di-n-propytamine	2-Nitrophenol	Hexachlorobutadiene
Hexachloroethane		Hexachlorocyclopentadiene
Nitrobenzene		4-Chloro-3-methylphenol
2,6-Dinitrotoluene		2,4,6-Trichlorophenol
2,4-Dinitrotoluene		2,4,5-Trichlorophenol
N-Nitrosodiphenylamine		1,2,4,5-Tetrachlorobenzene
		*Pentachlorophenol
		2,3,4,6-Tetrachlorophenol
Dimethylphthalate-d ₆ (DMC-10)	Acenaphthylene-da(DMC-11)	4-Nitrophenol-d ₄ (DMC-12)
Caprolactam	*Naphthalene	2-Nitroaniline
1,1'-Biphenyl	*2-Methylnaphthalene	3-Nitroaniline
Dimethylphthalate	2-Chloronaphthalene	2,4-Dinitrophenol
Diethylphthalate	*Acenaphthylene	4-Nitrophenol
Di-n-butylphthalate	*Acenaphthene	4-Nitroaniline
Butylbenzylphthalate	5.7	
Bis(2-ethylhexyl) phthalate		
Di-n-octylphthalate		

Fluorene-d ₁₀ (DMC-13)	4,6-Dinitro-2-methylphenol-d ₂ (DMC-14)	Anthracene-d ₁₀ (DMC-15)
Dibenzofuran *Fluorene 4-Chlorophenyl-phenylether 4-Bromophenyl-phenylether	4,6-Dinitro-2-methylphenol	Hexachlorobenzene Atrazine *Phenanthrene *Aothracene
Carbazole Pyrene-d ₁₀ (DMC-16)	Benzo(a)pyrene-d ₁₂ (DMC-17)	
*Fluoranthene *Pyrene *Benzo(a)anthracene *Chrysene	3,3'-Dichlorobenzidine *Benzo(b)fluoranthene *Benzo(k)fluoranthene *Benzo(a)pyrene *Indeno(1,2,3-cd)pyrene *Dibenzo(a,h)anthracene *Benzo(a,h,i)perylene	

^{*}Included in optional Target Analyte List (TAL) of PAHs and PCP only.

Table 9. Semivolatile SIM DMCs and the Associated Target Analytes

	., .
Fluoranthene-d10 (DMC-1)	2-Methylnaphthalene-d10 (DMC-2)
Fluoranthene	Naphthälene
Pyrene	2-Methylnaphthalene
Benzo(a)anthracene	Acenaphthylene
Chrysene	Acenaphthene
Benzo(b)fluoranthene	Fluorene
Benzo(k)fluoranthene	Pentachlorophenol
Benzo(a)pyrene	Phenanthrene
Indeno(1,2,3-ed)pyrene	Anthracene
Dibenzo(a,h)anthracene	
Benzo(g,h,i)perylene	

All criteria were metX	
Criteria were not met	
and/or see below	

VII. A MATRIX SPIKE/MATRIX SPIKE DUPLICATE (MS/MSD)

This data is generated to determine long term precision and accuracy in the analytical method for various matrices. This data alone cannot be used to evaluate the precision and accuracy of individual samples. If any % R in the MS or MSD falls outside the designated range, the reviewer should determine if there are matrix effects, i.e. LCS data are within the QC limits but MS/MSD data are outside QC limit.

1. MS/MSD Recoveries and Precision Criteria

The laboratory should use one MS and a duplicate analysis of an unspiked field sample if target analytes are expected in the sample. If target analytes are not expected, MS/MSD should be analyzed.

NOTES:

Data for MS and MSDs will not be present unless requested by the Region. Notify the Contract Laboratory COR if a field or trip blank was used for the MS and MSD.

For a Matrix Spike that does not meet criteria, apply the action to only the field sample used to prepare the Matrix Spike sample. If it is clearly stated in the data validation materials that the samples were taken through incremental sampling or some other method guaranteeing the homogeneity of the sample group, then the entire sample group may be qualified.

List the %Rs, RPD of the compounds which do not meet the criteria.

Sample ID:JC22506-9	Matrix/Level:Groundwater
Sample ID:JC22613-2	Matrix/Level:Groundwater
Sample ID:JC22506-8_(SIM)	Matrix/Level:Groundwater
Sample ID:JC22515-5_(SIM)	Matrix/Level:Groundwater

Note: MS/MSD % recoveries and RPD within laboratory control limits.

- * QC limits are laboratory in-house performance criteria, LL = lower limit, UL = upper limit.
- * If QC limits are not available, use limits of 70 130 %.

Actions:

1	QUALITY	%R < LL	%R > UL
	Positive results	J	J
	Nondetects results	R	Accept

MS/MSD criteria apply only to the unspiked sample, its dilutions, and the associated MS/MSD samples:

If the % R for the affected compounds were < LL (or 70 %), qualify positive results (J) and nondetects (UJ).

If the % R for the affected compounds were > UL (or 130 %), only qualify positive results (J). If 25 % or more of all MS/MSD %R were < LL (or 70 %) or if two or more MS/MSD %Rs were < 10%, qualify all positive results (J) and reject nondetects (R).

A separate worksheet should be used for each MS/MSD pair.

All criteria were met _X
Criteria were not met
and/or see below

INTERNAL STANDARD PERFORMANCE

The assessment of the internal standard (IS) parameter is used to assist the data reviewer in determining the condition of the analytical instrumentation.

List the internal standard area of samples which do not meet the criteria.

DATE SAMPLE ID IS OUT

IS AREA ACCEPTABLE RANGE

ACTION

Internal area meets the required criteria of batch samples corresponding to this data package.

Action:

- 1. If an internal standard area count for a sample or blank is greater than 200.0% of the area for the associated standard (opening CCV or mid-point standard from initial calibration) (see Table 10 below):
 - a. Qualify detects for compounds quantitated using that internal standard as estimated low (J-).
 - b. Do not qualify non-detected associated compounds.
- 2. If an internal standard area count for a sample or blank is less than 20.0% of the area for the associated standard (opening CCV or mid-point standard from initial calibration):
 - Qualify detects for compounds quantitated using that internal standard as estimated high (J+).
 - b. Qualify non-detected associated compounds as unusable (R).
- 3. If an internal standard area count for a sample or blank is greater than or equal to 50.0%, and less than or equal to 200% of the area for the associated standard opening CCV or mid-point standard from initial calibration, no qualification of the data is necessary.
- 4. If an internal standard RT varies by more than 10.0 seconds: Examine the chromatographic profile for that sample to determine if any false positives or negatives exist. For shifts of a large magnitude, the reviewer may consider partial or total rejection of the data for that sample fraction. Detects should not need to be qualified as unusable (R) if the mass spectral criteria are met.
- 5. If an internal standard RT varies by less than or equal to 10.0 seconds, no qualification of the data is necessary.

Note: Inform the Contract Laboratory Program Project Officer (CLP PO) if the internal standard performance criteria are grossly exceeded. Note in the Data Review Narrative potential effects on the data resulting from unacceptable internal standard performance.

State in the Data Review Narrative if the required internal standard compounds are not added to a sample or blank or if the required internal standard compound is not analyzed at the specified concentration.

Actions:

Table 10. Internal Standard Actions for Semivolatile Analysis

Criteria	Action		
Сгиена	Detect	Non-detect	
Area response < 20% of the opening CCV or mid-point standard CS3 from ICAL	11	R	
20% \le Area response \le 50% of the opening CCV or mid-point standard CS3 from ICAL	Jŧ	ÜJ	
50% < Area response < 200% of the opening CCV or mid-point standard CS3 from ICAL	No qualification	No qualification	
Area response > 200% of the opening CCV or mid-point standard CS3 from ICAL.	J-	No qualification	
RT shift between sample/blank and opening CCV or mid-point standard CS3 from ICAL > 10.0 seconds	R	R	
RT shift between sample/blank and opening CCV or mid-point standard CS3 from ICAL < 10.0 seconds	No qualification	No qualification	

		Alt criteria were metX Criteria were not met and/or see below
TARGET COM	MPOUND IDENTIFICATION	
Criteria:		
		oounds within ±0.06 RRT units of the standard CV) or mid-point standard from the initial Yes? or No?
List compound	ds not meeting the criteria described above:	
Sample ID	Compounds	Actions
spectrum from	n the associated calibration standard (openust match according to the following criteria All ions present in the standard mass spinust be present in the sample spectrum. The relative intensities of these ions mus sample spectra (e.g., for an ion with an the corresponding sample ion abundance lons present at greater than 10% in the	ectrum at a relative intensity greater than 10% tagree within ±20% between the standard and abundance of 50% in the standard spectrum,
List compound	ds not meeting the criteria described above:	
Sample ID	Compounds	Actions
ldentified_co	mpounds_meet_the_required_criteria	

Action:

- The application of qualitative criteria for GC/MS analysis of target compounds requires
 professional judgment. It is up to the reviewer's discretion to obtain additional information from
 the laboratory. If it is determined that incorrect identifications were made, qualify all such data
 as unusable (R).
- Use professional judgment to qualify the data if it is determined that cross-contamination has occurred.
- 3. Note in the Data Review Narrative any changes made to the reported compounds or concerns regarding target compound identifications. Note, for Contract Laboratory COR action, the necessity for numerous or significant changes.

TENTATIVELY IDENTIFIED COMPOUNDS (TICS)

NOTE: Tentatively identified compounds should only be evaluated when requested by a party from outside of the Hazardous Waste Support Section (HWSS).

	-	_
ict.	-14	/ 'c
131	- 1 - 3	1

Sample ID	Compound	Sample ID	Compound
=======================================			
	*******		- 10 Na

Action:

- 1. Qualify all TIC results for which there is presumptive evidence of a match (e.g. greater than or equal to 85% match) as tentatively identified (NJ), with approximated concentrations. TICs labeled "unknown" are qualified as estimated (J).
- 2. General actions related to the review of TIC results are as follows:
 - a. If it is determined that a tentative identification of a non-target compound is unacceptable, change the tentative identification to "unknown" or another appropriate identification, and qualify the result as estimated (J).
 - b. If all contractually-required peaks were not library searched and quantitated, the Region's designated representative may request these data from the laboratory.
- In deciding whether a library search result for a TIC represents a reasonable identification, use professional judgment. If there is more than one possible match, report the result as "either compound X or compound Y". If there is a lack of isomer specificity, change the TIC result to a nonspecific isomer result (e.g., 1,3,5-trimethyl benzene to trimethyl benzene isomer) or to a compound class (e.g., 2-methyl, 3-ethyl benzene to a substituted aromatic compound).
- 4. The reviewer may elect to report all similar compounds as a total (e.g., all alkanes may be summarized and reported as total hydrocarbons).

- 5. Target compounds from other fractions and suspected laboratory contaminants should be marked as "non-reportable".
- 6. Other Case factors may influence TIC judgments. If a sample TIC match is poor, but other samples have a TIC with a valid library match, similar RRT, and the same ions, infer identification information from the other sample TIC results.
- 7. Note in the Data Review Narrative any changes made to the reported data or any concerns regarding TIC identifications.
- 8. Note, for Contract Laboratory COR action, failure to properly evaluate and report TICs

All criteria were met_	_X	
Criteria were not met		
and/or see below		

SAMPLE QUANTITATION AND REPORTED CONTRACT REQUIRED QUANTITATION LIMITS (CRQLS)

Action:

- 1. When a sample is analyzed at more than one dilution, the lower CRQL are used unless a QC exceedance dictates the use of higher CRQLs from the diluted sample. Samples reported with an "E" qualifier should be reported from the diluted sample.
- 2. If any discrepancies are found, the Region's designated representative may contact the laboratory to obtain additional information that could resolve any differences. If a discrepancy remains unresolved, the reviewer must use professional judgment to decide which value is the most accurate. Under these circumstances, the reviewer may determine that qualification of data is warranted. Note in the Data Review Narrative a description of the reasons for data qualification and the qualification that is applied to the data.
- 3. For non-aqueous samples, if the solids is less than 10.0%, use professional judgment for both detects and non-detects. If the percent solid for a soil sample is greater than or equal to 10.0% and less than 30.0%, use professional judgment to qualify detects and non-detects. If the percent solid for a soil sample is greater than or equal to 30.0%, detects and non-detects should not be qualified (see Table 11).
- 4. Note, for Contract Laboratory COR action, numerous or significant failures to accurately quantify the target compounds or to properly evaluate and adjust CRQLs.
- 5. Results between MDL and CRQL should be qualified as estimated "J".
- 6. Results < MDL should be reported at the CRQL and qualified "U". MDLs themselves should not be reported.

Table 11. Percent Solids Actions for Semivolatile Analysis for Non-Aqueous Samples

Criteria	Action		
Criteria	Detects	Non-detects	
%Solids < 10.0%	Use professional judgment	Use professional judgment	
10.0% ≤ %Solids ≤ 30.0%	Use professional judgment	Use professional judgment	
%Solids > 30.0%	No qualification	No qualification	

SAMPLE QUANTITATION

The sample quantitation evaluation is to verify laboratory quantitation results. In the space below, please show a minimum of one sample calculation:

QUANTITATION LIMITS

A. Dilution performed

SAMPLE ID	DILUTION	REASON FOR DILUTION
A		
-		
, <u></u>	-	

				Crite	riteria were metX ria were not met or see below	
FIELD DUPLICATE	PRECIS	SION				
Sample IDs				Ма	trix:	
analyses measure laboratory duplicate will have a greater field duplicate samp The project QAPP s Suggested criteria:	both fieles which variance les. should be if large	d and lab preconly laboratory than water make reviewed for p RPD (> 50 %)	ision; therefore, the r r performance. It is al atrices due to difficult roject-specific informa	results may lso expected ies associat ation. identification	of overall precision. The have more variability of that soil duplicate reset with collecting identification of the samples and abled.	than sults tical
COMPOUND SQL SAMPLE DUPLICATE RPD ACTION ug/L CONC.						
					SD % recovery RPD ected target analytes	

All criteria were met _	_X	_
Criteria were not met		
and/or see below	_	

OTHER ISSUES

A.	System Per	formance	
List sa	amples qualifie	ed based on the degradation of system	performance during simple analysis:
Samp	le ID	Comments	Actions
Action	n:		
during	sample ana		nined that system performance has degraded y Program COR any action as a result of cted the data.
B.	Overall Asse	essment of Data	
List sa	amples qualifie	ed based on other issues:	
Samp	le ID	Comments	Actions
			_dataResults_are_valid_and_can_be_used n_below
Note:	JC22613-1	and JC22613-2: There are compound	s in BS were outside in house QC limits. The

Action:

1. Use professional judgment to determine if there is any need to qualify data which were not qualified based on the Quality Control (QC) criteria previously discussed.

results confirmed by re-extraction outside the holding time.

Write a brief narrative to give the user an indication of the analytical limitations of the data. Inform the Contract Laboratory COR the action, any inconsistency of the data with the Sample Delivery Group (SDG) Narrative. If sufficient information on the intended use and required quality of the data is available, the reviewer should include their assessment of the usability of the data within the given context. This may be used as part of a formal Data Quality Assessment (DQA).

- 3. Sometimes, due to dilutions, re-analysis or SIM/Scan runs are being performed, there will be multiple results for a single analyte from a single sample. The following criteria and professional judgment are used to determine which result should be reported:
 - The analysis with the lower CRQL
 - The analysis with the better QC results
 - The analysis with the higher results

EXECUTIVE NARRATIVE

SDG No:

JC22613

Laboratory:

Accutest, Florida

Analysis:

SW846-8015C

Number of Samples:

Я

Location:

BMSMC, Building 5 Area

Humacao, PR

SUMMARY:

Eight (8) samples were analyzed for the low molecular weight alcohols (LMWAs) list following method SW846-8015C. The sample results were assessed according to USEPA data validation guidance documents in the following order of precedence: "Test Methods for Evaluating Solid Waste, Physical/Chemical Methods SW-846 (Final Update III, December 1996)," specifically for Methods 8000/8015C are utilized. The QC criteria and data validation actions listed on the data review worksheets are from the primary guidance document, unless otherwise noted.

Results are valid and can be used for decision making purposes.

Critical issues:

None

Major:

None

Minor:

None

Critical findings:

None

Major findings:

None

Minor findings:

None

COMMENTS:

Results are valid and can be used for decision making purposes.

Reviewers Name:

Rafael Infante

Chemist License 1888

Signature:

July 21, 2016

Date:

SAMPLE ORGANIC DATA SAMPLE SUMMARY

Sample ID: JC22613-1

Sample location: BMSMC Building 5 Area

Sampling date: 6/16/2016

Matrix: Groundwater

METHOD: 8015C

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Ethanol	100	ug/l	1.0	-	U	Yes
Isobutyl Alcohol	100	ug/l	1.0	•	U	Yes
Isopropyl Alcohol	100	ug/l	1.0	-	Ü	Yes
n-Propyl Alcohol	100	ug/l	1.0	-	U	Yes
n-Butyl Alcohol	100	ug/l	1.0	-	U	Yes
sec-Butyl Alcohol	100	ug/l	1.0	-	U	Yes
Methanol	200	ug/l	1.0	-	IJ	Yes

Sample ID: JC22613-2

Sample location: BMSMC Building 5 Area

Sampling date: 6/17/2016

Matrix: Groundwater

METHOD: 8015C

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Ethanol	100	ug/l	1.0	-	U	Yes
Isobutyl Alcohol	100	ug/l	1.0	-	U	Yes
Isopropyl Alcohol	100	ug/l	1.0	-	U	Yes
n-Propyl Alcohol	100	ug/l	1.0	-	U	Yes
n-Butyl Alcohol	100	ug/l	1.0	-	U	Yes
sec-Butyl Alcohol	100	ug/l	1.0	-	U	Yes
Methanol	200	ug/l	1.0	-	U	Yes

Sample ID: JC22613-3

Sample location: BMSMC Building 5 Area

Sampling date: 6/17/2001 Matrix: Groundwater

METHOD: 8015C

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Ethanol	100	ug/l	1.0	•	U	Yes
Isobutyl Alcohol	100	ug/l	1.0	-	U	Yes
Isopropyl Alcohol	100	ug/l	1.0	-	U	Yes
n-Propyl Alcohol	100	ug/l	1.0	~	U	Yes
n-Butyl Alcohoi	100	ug/l	1.0	-	U	Yes
sec-Butyl Alcohol	100	ug/l	1.0	-	U	Yes
Methanol	200	ug/l	1.0	-	U	Yes

Sample ID: JC22613-4

Sample location: BMSMC Building 5 Area

Sampling date: 6/17/2016 Matrix: Groundwater

METHOD: 8015C

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Ethanol	100	ug/l	1.0	-	U	Yes
Isobutyl Alcohol	100	ug/l	1.0	•	U	Yes
Isopropyl Alcohol	100	ug/l	1.0	-	U	Yes
n-Propyl Alcohol	100	ug/l	1.0	-	U	Yes
n-Butyl Alcohol	100	ug/l	1.0	•	U	Yes
sec-Butyl Alcohol	100	ug/l	1.0	-	บ	Yes
Methanol	200	ug/l	1.0	-	U	Yes

Sample ID: JC22613-5

Sample location: BMSMC Building 5 Area

Sampling date: 6/17/2016 Matrix: Groundwater

METHOD: 8015C

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Ethanol	100	ug/l	1.0	-	U	Yes
isobutyl Alcohol	100	ug/l	1.0	-	U	Yes
isopropyl Alcohol	100	ug/l	1.0	-	U	Yes
n-Propyl Alcohol	100	ug/l	1.0	-	U	Yes
n-Butyl Alcohol	100	ug/l	1.0	-	U	Yes
sec-Butyl Alcohol	100	ug/l	1.0	-	U	Yes
Methanol	200	ug/l	1.0	-	U	Yes

Sample ID: JC22613-6

Sample location: BMSMC Building 5 Area

Sampling date: 6/17/2016 Matrix: Groundwater

METHOD: 8015C

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Ethanol	100	ug/l	1.0	-	U	Yes
Isobutyl Alcohol	100	ug/l	1.0	-	U	Yes
Isopropyl Alcohol	100	ug/l	1.0	-	U	Yes
n-Propyl Alcohol	100	ug/l	1.0	-	U	Yes
n-Butyl Alcohol	100	ug/l	1.0	•	U	Yes
sec-Butyl Alcohol	100	ug/l	1.0	-	υ	Yes
Methanol	200	ug/l	1.0	-	υ	Yes

Sample ID: JC22613-7

Sample location: BMSMC Building 5 Area

Sampling date: 6/20/2016 Matrix: Groundwater

METHOD: 8015C

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Ethanol	100	ug/l	1.0	-	υ	Yes
Isobutyl Alcohol	100	ug/l	1.0	-	U	Yes
Isopropyi Alcohol	100	ug/l	1.0	-	U	Yes
n-Propyl Alcohol	100	ug/l	1.0	-	U	Yes
n-Butyl Alcohol	100	ug/l	1.0	-	U	Yes
sec-Butyl Alcohol	100	ug/l	1.0	-	U	Yes
Methanol	200	ug/l	1.0	-	U	Yes

Sample ID: JC22613-8

Sample location: BMSMC Building 5 Area

Sampling date: 6/20/2016 Matrix: Groundwater

METHOD: 8015C

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Ethanol	100	ug/l	1.0	•	U	Yes
Isobutyl Alcohol	100	ug/l	1.0		U	Yes
Isopropyl Alcohol	100	ug/l	1.0	7.	U	Yes
n-Propyl Alcohol	100	ug/l	1.0		U	Yes
n-Butyl Alcohol	100	ug/l	1.0		U	Yes
sec-Butyl Alcohol	100	ug/l	1.0	-	U	Yes
Methanol	200	ug/l	1.0		U	Yes

	Project Number:JC22613
	Date:06/16-20/2016
	Shipping Date:06/20/2016
	EPA Region:2_
REVIEW OF VOLATILE (The following guidelines for evaluating volatile organics were document will assist the reviewer in using professional judgerving the needs of the data users. The sample results guidance documents in the following order of preceded Physical/Chemical Methods SW-846 (Final Update III, Decentialized. The QC criteria and data validation actions listed guidance document, unless otherwise noted. The hardcopied (laboratory name) _Accutest	ORGANIC PACKAGE re created to delineate required validation actions. This dgment to make more informed decision and in better were assessed according to USEPA data validation lence: "Test Methods for Evaluating Solid Waste, ember 1996)," specifically for Methods 8000/8015C are d on the data review worksheets are from the primary data package received has been reviewed The modified data review for VOCs included:
No. of Samples:8	
Frield blank No.: Equipment blank No.: X Data Completeness X Holding Times N/A_ GC/MS Tuning N/A_ Internal Standard Performance X Blanks X Blanks X Surrogate Recoveries X Matrix Spike/Matrix Spike Duplicate Overall Comments: _Low_molecular_weight_alcompany contents and comments and comments are contents and conte	X Laboratory Control SpikesX Field DuplicatesX CalibrationsX Compound IdentificationsX Compound QuantitationX Quantitation Limits
Overali Comments:_Low_molecular_weight_aid	:onois_by_5vv-846_8015C
Definition of Qualifiers: J- Estimated results J- Compound not detected R- Rejected data JJ- Estimated nondetect Reviewer: Date: July 21, 2016	

DATA COMPLETENESS

MISSING INFORMATION	DATE LAB. CONTACTED	DATE RECEIVED
		· · · · · · · · · · · · · · · · · · ·
	4. CD N 197.350	
-		
	>	
* The Control of the		
	7	
		(19
	4	
	-	

All criteria were met _	_X
Criteria were not met	
and/or see below	

HOLDING TIMES

The objective of this parameter is to ascertain the validity of the results based on the holding time of the sample from time of collection to the time of analysis.

Complete table for all samples and note the analysis and/or preservation not within criteria

DATE SAMPLED	DATE ANALYZED	рН	ACTION
lyzed within the red	commended method ho	olding ti	me. All samples properly
			DATE SAMPLED DATE ANALYZED PH lyzed within the recommended method holding ti

Criteria

Aqueous samples – 14 days from sample collection for preserved samples (pH \leq 2, 4°C), no air bubbles. Aqueous samples – 7 days from sample collection for unpreserved samples, 4°C, no air bubbles.

Soil samples- 7 days from sample collection.

Cooler temperature (Criteria: 4 ± 2 °C): 3.9°C

Actions

If the VOCs vial(s) have air bubbles, estimate positive results (J) and reject nondetects (R).

If the % solids of soil samples is 10-50%, estimates positive results (J) and nondetects (UJ)

If the % solid of soil samples is < 10%, estimate positive results (J) and reject nondetects (R).

If holding times are exceeded but < 14 days beyond criteria, estimate positive results (J) and nondetects (UJ).

If holding times are exceeded but < 28 days beyond criteria, estimate positive results (J) and reject nondetects (R).

If holding times are grossly exceeded (> 28 days beyond criteria), reject all results (R).

If samples were not iced or if the ice were melted (> 10°C), estimate positive results (J) and nondetects (UJ).

	All criteria we	re met	N/A_
Criteria	were not met s	ee below	

GC/MS TUNING

The assessment of the tuning results is to determine if the sample instrumentation is within the standard uning QC limits
N/A The BFB performance results were reviewed and found to be within the specified criteria.
N/A BFB tuning was performed for every 12 hours of sample analysis.
f no, use professional judgment to determine whether the associated data should be accepted, qualified or rejected.
ist the samples affected:

If mass calibration is in error, all associated data are rejected.

All criteria were met _X_	
Criteria were not met	
and/or see below	

CALIBRATION VERIFICATION

Compliance requirements for satisfactory instrument calibration are established to ensure that the instrument is capable of producing and maintaining acceptable quantitative data.

	Date	e of initial calibration:	05/17/16		
	Date	es of continuing calibra	tion:05/17/16	(initial);_06/23/16	
	Date	es of final calibration ve	erification:06/23/16_		
	Instrument ID number: GCGH				
	Mat	rix/Level:	Aqueous/low		
			•		
DATE	LAB FILE ID#	CRITERIA OUT	COMPOUND	SAMPLES	
		RFs, %RSD, %D, r		AFFECTED	

	DATE	LAB FILE ID#	CRITERIA OUT RFs, %RSD, %D, r	COMPOUND	SAMPLES AFFECTED
-					C
ŀ					

Note: Initial, continuing, and final calibration verifications meets method specific criteria.

Criteria

г

All RFs must be > 0.05 regardless of method requirements for SPCC.

All %RSD must be < 15 % regardless of method requirements for CCC.

All %Ds must be \leq 20% regardless of method requirements for CCC.

It should be noted that Region 2 SOP HW-24 does not specify criterion for the curve correlation coefficient (r). A limit for r of \geq 0.995 has therefore been utilized as professional judgment.

Actions

If any compound has an initial RF or a continuing RF of < 0.05, estimate positive results (J) and reject nondetects (R), regardless of method requirements.

If any compound has a %RSD > 15%, estimate positive results (J) and use professional judgment to qualify nondetects.

If any compound has a %RSD > 90%, estimate positive results (J) and reject nondetects (R).

If any compound has a % D > 20%, estimate positive results (J) and reject nondetects (R).

If any compound has a % D > 20%, estimate positive results (J) and nondetects (UJ).

If any compound has a % D > 90%, estimate positive results (J) and reject nondetects (R).

If any compound has r < 0.995, estimate positive results and nondetects.

A separate worksheet should be filled for each initial curve

All criteria were met _X
Criteria were not met
and/or see below

V A. BLANK ANALYSIS RESULTS (Sections 1 & 2)

The assessment of the blank analysis results is to determine the existence and magnitude of contamination problems. The criteria for evaluation of blanks apply only to blanks associated with the samples, including trip, equipment, and laboratory blanks. If problems with any blanks exist, all data associated with the case must be carefully evaluated to determine whether or not there is an inherent variability in the data for the case, or if the problem is an isolated occurrence not affecting other data.

List the contamination in the blanks below. High and low levels blanks must be treated separately.

Laboratory blanks

DATE ANALYZED	LAB ID	LEVEL/ MATRIX	COMPOUND	CONCENTRATION UNITS
All_method			fic_criteria	
Field/Equipment		LEVEL	COMPOUND	CONCENTRATION
ANALYZED _No_field/trip/ed	uipment_blank	MATRIX s_included_in_	this_data_package	UNITS
		1		

All criteria were met _X
Critena were not met
and/or see below

VB. BLANK ANALYSIS RESULTS (Section 3)

Blank Actions

Action Levels (ALs) should be based upon the highest concentration of contaminant determined in any blank. Do not qualify any blank with another blank. The ALs for samples which have been diluted should be corrected for the sample dilution factor and/or % moisture, where applicable. No positive sample results should be reported unless the concentration of the compound in the samples exceeds the ALs:

ALs = 10x the amount of common contaminants (methylene chloride, acetone, 2-butanone, and toluene) ALs = 5x for any other compounds

Specific actions are as follows:

If the concentration is < sample quantitation limit (SQL) and \le AL, report the compound as not detected (U) at the SQL.

If the concentration is \geq SQL but \leq AL, report the compound as not detected (U) at the reported concentration.

If the concentration is \geq SQL and > AL, report the concentration unqualified.

Notes:

High and low level blanks must be treated separately Compounds qualified "U" for blank contamination are still considered "hits" when qualifying for calibration

Compounds qualified "O" for blank contamination are still considered "hits" when qualifying for calibration criteria.

CONTAMINATION SOURCE/LEVEL	COMPOUND	CONC/UNITS	AL/UNITS	SQL	AFFECTED SAMPLES
			ļ		Alder-
			and a		
		4			
100					
E					<u> </u>

All criteria were metX	
Criteria were not met	
and/or see below	

SURROGATE SPIKE RECOVERIES

Laboratory performance of individual samples is established by evaluation of surrogate spike recoveries. All samples are spiked with surrogate compounds prior to sample analysis. The accuracy of the analysis is measured by the surrogate percent recovery. Since the effects of the sample matrix are frequently outside the control of the laboratory and may present relatively unique problems, the validation of data is frequently subjective and demands analytical experience and professional judgment.

List the percent recoveries (%Rs) which do not meet the criteria for surrogate recovery. Matrix: solid/aqueous

SAMPLE ID		SURROGATE COMPOUND A			ACTION
	Hexanol	DBFM	TOL-d8	BFB	
_All_surrogate_	recoveries_withi	n_laboratory_co	ntrol_limits		
	01 12 10 16				
	JL56_to_	145to	to	to	
QC Limits* (Soli LL_to_U QC Limits* (Soli	Lto	to	to	to	_
LL_to_U	Lto	to	to	to	
	Dichloromethane nofluoromethane	-d4		3 = Toluene-d8 Bromofluorobenze	ne
* If QC lir				= lower limit, UL = aqueous and 70 -	
samples.					
Actions:					
QUALI	TY	%R < 10%	%R = 1	0% - LL %R >	· UL
Positive	e results	J	J	J	
Nonde	tects results	R	UJ	Acce	pt

Surrogate action should be applied:

If one or more surrogate in the VOC fraction is out of specification, but has a recovery of > 10%. If any one surrogate in a fraction shows < 10 % recovery.

All criteria were metX
Criteria were not met
and/or see below

VII. A MATRIX SPIKE/MATRIX SPIKE DUPLICATE (MS/MSD)

This data is generated to determine long term precision and accuracy in the analytical method for various matrices. This data alone cannot be used to evaluate the precision and accuracy of individual samples. If any % R in the MS or MSD falls outside the designated range, the reviewer should determine if there are matrix effects, i.e. LCS data are within the QC limits but MS/MSD data are outside QC limit.

1. MS/MSD Recoveries and Precision Criteria

The laboratory should use one MS and a duplicate analysis of an unspiked field sample if target analytes are expected in the sample. If target analytes are not expected, MS/MSD should be analyzed.

List the %Rs, RPD of the compounds which do not meet the criteria.

Sample ID:JC22613-1MS/-MSD				Matrix/Level:Groundwater/low	
MS OR MSD	COMPOUND	% R	RPD	QC LIMITS	ACTION
MS/MSD%_re	ecoveries_and_RPD_w	ithin_labo	oratory_	control_limits_ex	cept_for_the_following:
_MSD	lsobutyl_alcohol	134_%		69131	No_action

Note: No action, professional judgment. Analyte not detected in sample batch, non-detects are accepted.

- * QC limits are laboratory in-house performance criteria, LL = lower limit, UL = upper limit.
- * If QC limits are not available, use limits of 70 130 %.

Actions:

QUALITY	%R < LL	%R > UL
Positive results	J	J
Nondetects results	R	Accept

All criteria were met _X
Criteria were not met
and/or see below

MS/MSD criteria apply only to the unspiked sample, its dilutions, and the associated MS/MSD samples:

If the % R for the affected compounds were < LL (or 70 %), qualify positive results (J) and nondetects (UJ).

If the % R for the affected compounds were > UL (or 130 %), only qualify positive results (J). If 25 % or more of all MS/MSD %R were < LL (or 70 %) or if two or more MS/MSD %Rs were < 10%, qualify all positive results (J) and reject nondetects (R).

VII. B MATRIX SPIKE/MATRIX SPIKE DUPLICATE

MS/MSD - Unspiked Compounds

It should be noted that Region 2 SOP HW-24 does not specify a MS/MSD criteria for the unspiked compounds in the sample. A %RSD of < 50% has therefore been utilized as professional judgment.

If all target analytes were spiked in the MS/MSD, this review element is not applicable.

List the %RSD of the compounds which do not meet the criteria.

Sample ID:			Matrix/Le	vel/Unit	-	
COMPOUND	SAMPLE CONC.	MS CONC.	MSD CONC.	% RSD	ACTION	
					La Contraction of the Contractio	
		20,7 1.1 70				POST
	974			10 1 10 10 10 10 10 10 10 10 10 10 10 10		
<u>. </u>		The state of the s	······································			
priting.				***	· · · · · · · · · · · · · · · · · · ·	

Actions:

A separate worksheet should be used for each MS/MSD pair.

^{*} If the % RSD > 50, qualify the positive result in the unspiked samples as estimated (J).

^{*} If the % RSD is not calculated (NC) due to nondetected value, use professional judgment to qualify the data.

All criteria were met _X	
Criteria were not met	
and/or see below	

VIII. LABORATORY CONTROL SAMPLE (LCS) ANALYSIS

This data is generated to determine accuracy of the analytical method for various matrices.

1. LCS Recoveries Criteria

Where LCS spiked with the same analyte at the same concentrations as the MS/MSD? Yes or No. If no make note in data review memo.

List the %R of compounds which do not meet the criteria

	LCS ID	COMPOUND	% R	QC LIMIT	
Recoverie	s_within_labor	ratory_control_limits			
	7/23 .1-		W 30 12 13 13 13 13 13 13 13 13 13 13 13 13 13		-0.20
					_

- * QC limits are laboratory in-house performance criteria, LL = lower limit, UL = upper limit
- * If QC limits are not available, use limits of 70 130 %.

Actions:

QUALITY	%R < LL	%R > UL
Positive results	J	J
Nondetects results	R	Accept

All analytes in the associated sample results are qualified for the following criteria.

If 25 % of the LCS recoveries were < LL (or 70 %), qualify all positive results (j) and reject nondetects (R).

If two or more LCS were below 10 %, qualify all positive results as (J) and reject nondetects (R).

2. Frequency Criteria:

Where LCS analyzed at the required frequency and for each matrix? <u>Yes</u> or No. If no, the data may be affected. Use professional judgment to determine the severity of the effect and qualify data accordingly. Discuss any actions below and list the samples affected.

		All criteria were met Criteria were not met and/or see belowN/A
IX.	FIELD/LABORATORY DUPLICATE PRECISION	
	Sample IDs:	Matrix:

Field/laboratory duplicates samples may be taken and analyzed as an indication of overall precision. These analyses measure both field and lab precision; therefore, the results may have more variability than laboratory duplicates which only laboratory performance. It is also expected that soil duplicate results will have a greater variance than water matrices due to difficulties associated with collecting identical field duplicate samples.

The project QAPP should be reviewed for project-specific information. Suggested criteria: RPD \pm 30% for aqueous samples, RPD \pm 50 % for solid samples. If both samples and duplicate are <5 SQL, the RPD criteria is doubled.

COMPOUND	SQL	SAMPLE CONC.	DUPLICATE CONC.	RPD	ACTION	
No field/laboratory duplicate analyzed with this data package. MS/MSD % recovery RPD used to assess precision. RPD within laboratory, generally acceptable and guidance document performance criteria control limits.						

Actions:

Qualify as estimated positive results (J) and nondetects (UJ) for the compound that exceeded the above criteria. For organics, only the sample and duplicate will be qualified.

If an RPD cannot be calculated because one or both of the sample results is not detected, the following actions apply:

If one sample result is not detected and the other is greater than 5x the SQL qualify (J/UJ).

If one sample value is not detected and the other is greater than 5x the SQL and the SQLs for the sample and duplicate are significantly different, use professional judgment to determine if qualification is appropriate.

If one sample value is not detected and the other is less than 5x, use professional judgment to determine if qualification is appropriate.

If both sample and duplicate results are not detected, no action is needed.

All criteria were met __N/A__ Criteria were not met and/or see below ____

X. INTERNAL STANDARD PERFORMANCE

The assessment of the internal standard (IS) parameter is used to assist the data reviewer in determining the condition of the analytical instrumentation.

List the internal standard area of samples which do not meet the criteria.

- * Area of +100% or -50% of the IS area in the associated calibration standard.
- * Retention time (RT) within 30 seconds of the IS area in the associated calibration standard.

		RANGE	
	3.250 <u>1. 3.</u>		- =
<u>.</u>			
 		SEE SEE	
	E STATE OF THE STA		

Actions:

1. IS actions should be applied to the compound quantitated with the out-of-control ISs

QUALITY	IS AREA < -25%	IS AREA = -25 %	IS AREA > + 100%
		TO – 50%	
Positive results	J	J	J
Nondetected results	R	UJ	ACCEPT

 If a IS retention time varies more than 30 seconds, the chromatographic profile for that sample must be examined to determine if any false positive or negative exists. For shifts of a large magnitude, the reviewer may consider partial or total rejection of the data for the sample fraction.

All criteria were metX_	
Criteria were not met	
and/or see below	

XII. SAMPLE QUANTITATION

The sample quantitation evaluation is to verify laboratory quantitation results. In the space below, please show a minimum of one sample calculation:

JC22613-1

Hexanol

RF = 67.60

[] = (264452)/(67.60)

= 3,912 ppm OK

All criteria were metX
Criteria were not met
and/or see below

XII. QUANTITATION LIMITS

A. Dilution performed

SAMPLE ID	DILUTION FACTOR	REASON FOR DILUTION		
		- 01		

B.	Percent	Solids	

List samples wh	ich have ≤ {	50 % solid	ds			
					 - 688	- 10
				40000		

Actions:

If the % solids of a soil sample is 10-50%, estimate positive results (J) and nondetects (UJ)

If the % solids of a soil sample is < 10%, estimate positive results (J) and reject nondetects (R) $\,$

EXECUTIVE NARRATIVE

SDG No:

JC22613

Laboratory:

Accutest, New Jersey

Analysis:

SW846-8081B

Number of Samples:

2

Location:

BMSMC, Building 5 Area

Humacao, PR

SUMMARY:

Two (2) samples were analyzed for selected pesticides following method SW846-8081B. The sample results were assessed according to USEPA data validation guidance documents in the following order of precedence *Hazardous Waste Support Section SOP No. HW-36A, Revision 0, June, 2015. SOM02.2. Pesticide Data Validation.* The QC criteria and data validation actions listed on the data review worksheets are from the primary guidance document, unless otherwise noted.

Results are valid and can be used for decision making purposes.

Critical issues:

None

Major:

None

Minor:

None

Critical findings:

None

Major findings:

None

Minor findings:

1. Initial and initial calibration verification within the guidance document performance criteria. Continuing calibration % differences meet the performance criteria in at least one of the two columns. Final calibration verification not included in data package. No action taken, professional judgment.

carron, protossional judginoni

COMMENTS:

Results are valid and can be used for decision making purposes.

Reviewers Name:

Rafael Infante

Chemist License 1888

Signature:

Date:

July 21, 2016

SAMPLE ORGANIC DATA SAMPLE SUMMARY

....

Sample ID: JC22439-7

Sample location: BMSMC Building 5 Area

Sampling date: 20-Jun-16 Matrix: Groundwater

METHOD: 8081B

171211100	,					
Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Aldrin	0.010	ug/l	1	-	ប	Yes
alpha-BHC	0.010	ug/l	1	-	U	Yes
beta-BHC	0.010	ug/l	1	-//	U	Yes
delta-BHC	0.010	ug/l	1	-	U	Yes
gamma-BHC (Lindane)	0.010	ug/l	1	-	U	Yes
alpha-Chlordane	0.010	ug/l	1	-	U	Yes
gamma-Chlordane	0.010	ug/l	1	-	U	Yes
Dieldrin	0.010	ug/l	1	-	U	Yes
4,4'-DDD	0.010	ug/l	1	-	U	Yes
4,4'-DDE	0.010	ug/l	1	-	U	Yes
4,4'-DDT	0.010	ug/l	1	-	U	Yes
Endrin	0.010	ug/l	1	-	U	Yes
Endosulfan sulfate	0.010	ug/l	1	-	U	Yes
Endrin aldehyde	0.010	ug/l	1	-	U	Yes
Endrin ketone	0.010	ug/l	1	-	U	Yes
Endosulfan-l	0.010	ug/l	1	-	ប	Yes
Endosulfan-II	0.010	ug/l	1	-	U	Yes
Heptachlor	0.010	ug/l	1	-	U	Yes
Heptachlor epoxide	0.010	ug/l	1	-	U	Yes
Methoxychlor	0.020	ug/l	1	-	U	Yes
Toxaphene	0.25	ug/l	1	-	U	Yes

Sample ID: JC22439-2

....

Sample location: BMSMC Building 5 Area

Sampling date: 20-Jun-16 Matrix: Groundwater

METHOD: 8081B

****	100. 00020					
Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Aldrin	0.010	ug/l	1	-	U	Yes
alpha-BHC	0.010	ug/l	1	-	U	Yes
beta-BHC	0.010	ug/l	1	-	U	Yes
delta-BHC	0.010	ug/l	1	-	U	Yes
gamma-BHC (Lindane)	0.010	ug/l	1	_	IJ	Yes
alpha-Chlordane	0.010	ug/l	1	-	บ	Yes
gamma-Chlordane	0.010	ug/l	1	-	U	Yes
Dieldrin	0.010	ug/l	1	-	U	Yes
4,4'-DDD	0.010	ug/l	1	-	U	Yes
4,4'-DDE	0.010	ug/l	1	-	U	Yes
4,4'-DDT	0.010	ug/l	1	-	U	Yes
Endrin	0.010	ug/l	1	-	U	Yes
Endosulfan sulfate	0.010	ug/l	1	-	U	Yes
Endrin aldehyde	0.010	ug/l	1	-	IJ	Yes
Endrin ketone	0.010	ug/l	1	-	U	Yes
Endosulfan-i	0.010	ug/l	1	-	IJ	Yes
Endosulfan-II	0.010	ug/l	1	-	U	Yes
Heptachlor	0.010	ug/l	1	-	υ	Yes
Heptachlor epoxide	0.010	ug/l	1	-	U	Yes
Methoxychlor	0.020	ug/l	1	-	U	Yes
Toxaphene	0.26	ug/l	1	-	U	Yes

	Project/Case Number:JC22613 Sampling Date:June_16-20,_2016 Shipping Date:June_20,_2016 EPA Region No.:2
REVIEW OF PESTICIDE ORGA	ANIC PACKAGE
The following guidelines for evaluating volatile required validation actions. This document will ass judgment to make more informed decision and in users. The sample results were assessed according documents in the following order of precedence Ha HW-36A, Revision 0, June, 2015. SOM02.2. Pesticided data validation actions listed on the data review guidance document, unless otherwise noted.	bist the reviewer in using professional better serving the needs of the data g to USEPA data validation guidance zardous Waste Support Section SOP No. Data Validation. The QC criteria and
The hardcopied (laboratory name) _Accutest	data package received has been rized. The data review for VOCs included:
Lab. Project/SDG No.:JC22613 No. of Samples:2	Sample matrix:Groundwater
Trip blank No.:	
X Data CompletenessX Holding TimesN/A GC/MS TuningX_ Internal Standard PerformanceX BlanksX Surrogate RecoveriesX Matrix Spike/Matrix Spike Duplicate	XLaboratory Control SpikesX Field DuplicatesX CalibrationsX Compound IdentificationsX Compound QuantitationX Quantitation Limits
Overall Comments:TCL_pesticides_list_by_SW846-80	81B
Definition of Qualifiers: J- Estimated results U- Compound not detected R- Rejected data UJ- Estimated nordeted Reviewer: Date: July 21, 2016	

DATA COMPLETENESS

MISSING INFORMATION	DATE LAB. CONTACTED	DATE RECEIVED
	7/4	114
-		
		···
		**
	-	
<u> </u>		
		7
		<u> </u>
		- "
	2.77	
	West and control	Δ.
194		

All criteria were metX
Criteria were not met
and/or see below

HOLDING TIMES

The objective of this parameter is to ascertain the validity of the results based on the holding time of the sample from time of collection to the time of analysis.

Complete table for all samples and note the analysis and/or preservation not within criteria

SAMPLE ID	DATE SAMPLED	DATE EXTRACTED/ANALYZED	ACTION
Samples properly pro	eserved.		

Preservatives:	Ali_samples_	_extracted_	and_	_analyzed_	_within_	the_	_required_	_criteria	

Criteria

Aqueous samples - seven (7) days from sample collection for extraction; 40 days from sample collection for analysis.

Non-aqueous samples – fourteen (14) days from sample collection for extraction; 40 days from sample collection for analysis.

Cooler temperature (Criteria: 4 ± 2 °C): 3.9°C - OK

Actions

Qualify aqueous sample results using preservation and technical holding time information as follows:

- a. If there is no evidence that the samples were properly preserved (T = 4° C \pm 2° C), and the samples were extracted or analyzed within the technical holding times, qualify detects as estimated (J) and non-detects as estimated (UJ).
- b. If there is no evidence that the samples were properly preserved (T = 4° C \pm 2° C), and the samples were extracted or analyzed outside the technical holding times, qualify detects as estimated (J) and non-detects as estimated (UJ).
- c. If the samples were properly preserved, and were extracted and analyzed within the technical holding times, no qualification of the data is necessary.
- d. If the samples were properly preserved, and were extracted or analyzed outside the technical holding times, qualify detects as estimated (J) and non-detects as estimated (UJ). Note in the Data Review Narrative that holding times were exceeded and the effect of exceeding the holding time on the resulting data.

- e. Use professional judgment to qualify samples whose temperature upon receipt at the laboratory is either below 2 degrees centigrade or above 6 degrees centigrade.
- f. If technical holding times are grossly exceeded, use professional judgment to qualify the data.

Qualify non-aqueous sample results using preservation and technical holding time information as follows:

- a. If there is no evidence that the samples were properly preserved (T = 4° C \pm 2° C), and the samples were extracted or analyzed within the technical holding time, qualify detects as estimated (J) and non-detects as estimated (UJ).
- b. If there is no evidence that the samples were properly preserved (T = 4° C \pm 2° C), and the samples were extracted or analyzed outside the technical holding time, qualify detects as estimated (J) and non-detects as estimated (UJ).
- c. If the samples were properly preserved, and were extracted and analyzed within the technical holding time, no qualification of the data is necessary.
- d. If the samples were properly preserved, and were extracted or analyzed outside the technical holding time, qualify detects as estimated (J) and non-detects as estimated (UJ). Note in the Data Review Narrative that holding times were exceeded and the effect of exceeding the holding time on the resulting data.
- e. Use professional judgment to qualify samples whose temperature upon receipt at the laboratory is either below 2 degrees centigrade or above 6 degrees centigrade.
- f. If technical holding times are grossly exceeded, use professional judgment to qualify the data.

	All criteria were met	X
Criteria	were not met see below	

GAS CHROMATOGRAPH WITH ELECTRON CAPTURE DETECTOR (GC/ECD) INSTRUMENT PERFORMANCE CHECK (SECTIONS 1 TO 5)

1. Resolution Check Mixture

Criteria

Is the resolution between two adjacent peaks in the Resolution Check Mixture C greater than or equal to 80.0% for all analytes for the primary column and greater than or equal to 50.0% for the confirmation column?

Yes? or No?

Is the resolution between two adjacent peaks in the Resolution Check Mixture (A and B) greater than or equal to 60.0%?

Yes? or No?

Note:

If resolution criteria are not met, the quantitative results may not be accurate due to inadequate resolution. Qualitative identifications may also be questionable if coelution exists.

Action

- a. Qualify detects for target compounds that were not adequately resolved as tentatively identified (NJ).
- b. Qualify non-detected compounds as unusable (R).

2. Performance Evaluation Mixture (PEM) Resolution Criteria

Criteria

Is PEM analysis performed at the required frequency (at the end of each pesticide initial calibration sequence and every 12 hours)?

Yes? or No?

Action

a. If PEM is not performed at the required frequency, qualify all associated sample and blank results as unusable (R).

Criteria

Is PEM % Resolution < 90%?

Yes? or No?

Action

- a. a. Qualify detects for target compounds that were not adequately resolved as tentatively identified (NJ).
- b. Qualify non-detected compounds as unusable (R).

	All criteria were met	_X
Criteria	were not met see below	7.11.30

3. PEM 4,4'-DDT Breakdown

Criteria

Is the PEM 4,4'-DDT % Breakdown >20.0% and 4,4'-DDT is detected?

Yes? or No?

Action

a. Qualify detects for 4,4'-DDT; detects for 4,4'-DDD; and detects for 4,4'-DDE as estimated (J)

Criteria

Is the PEM 4,4'-DDT % Breakdown >20.0% and 4,4'-DDT is not detected

Yes? or No?

Action

- a. Qualify non-detects for 4,4'- DDT as unusable (R)
- b. Qualify detects for 4.4'-DDD as tentatively identified (NJ)
- c. Qualify detects for 4,4'-DDE as tentatively identified (NJ)

4. PEM Endrin Breakdown

Criteria

Is the PEM Endrin % Breakdown >20.0% and Endrin is detected?

Yes? or No?

Action

a. Qualify detects for Endrin; detects for Endrin aldehyde; and detects for Endrin ketone as estimated (J)

Criteria

Is the PEM Endrin % Breakdown >20.0% and Endrin is not detected

Yes? or No?

Action

- a. Qualify non-detects for Endrin as unusable (R)
- b. Qualify detects for Endrin aldehyde as tentatively identified (NJ)
- c. Qualify detects for Endrin ketone as tentatively identified (NJ)

	All criteria were metX	
Criteria	were not met see below	

5. Mid-point Individual Standard Mixture Resolution -

Criteria

Is the resolution between two adjacent peaks in the Resolution Check Mixture C greater than or equal to 80.0% for all analytes for the primary column and greater than or equal to 50.0% for the confirmation column?

Yes? or No?

Is the resolution between two adjacent peaks in the Resolution Check Mixture (A and B) greater than or equal to 90.0%?

Yes? or No?

Note: If resolution criteria are not met, the quantitative results may not be accurate due to inadequate resolution. Qualitative identifications may also be questionable if coelution exists.

Action

- a. Qualify detects for target compounds that were not adequately resolved as tentatively identified (NJ).
- b. Qualify non-detected compounds as unusable (R).

Criteria

Is mid-point individual standard mixture analysis performed at the required frequency (every 12 hours)?

Yes? or No?

Action

a. If the mid-point individual standard mixture analysis is not performed at the required frequency, qualify all associated sample and blank results as unusable (R).

All criteria were metX
Criteria were not met
and/or see below

CALIBRATION VERIFICATION

Compliance requirements for satisfactory instrument calibration are established to ensure that the instrument is capable of producing and maintaining acceptable quantitative data.

Date of initial calibration:	06/10/16
Dates of initial calibration verific	ation:06/10/16
Dates of continuing calibration:_	06/27/16
Dates of final calibration	06/28/16
Instrument ID numbers:	_GC6G
Matrix/Level:	Aqueous/low

DATE	LAB	FILE	CRITERIA OUT	COMPOUND	SAMPLES AFFECTED
	ID#		RFs, %RSD, %D, r		
					ment performance criteria.
					ia in at least one of the two
columns	. Final c	alibratio			No action taken, professional
			าเกีย	ment.	
	_		Jaag	1110110	
			Judg		

Criteria

Are a five point calibration curve delivered with concentration levels as shown in Table 3 of SOP HW-36A, Revision 0, June, 2015?

Yes? or No?

Actions

If the standard concentrations listed in Table 3 are not used, use professional judgment to evaluate the effect on the data

Criteria

Are RT Windows calculated correctly?

Yes? or No?

All criteria were metX
Criteria were not met
and/or see below

Action

Recalculate the windows and use the corrected values for all evaluations.

Criteria

Are the Percent Relative Standard Deviation (%RSD) of the CFs for each of the single component target compounds less than or equal to 20.0%, except for alpha-BHC and delta-BHC?

Yes? or No?

Are the %RSD of the CFs for alpha-BHC and delta-BHC less than or equal to 25.0%. Yes? or No?

Is the %RSD of the CFs for each of the Toxaphene peaks must be < 30% when 5-point ICAL is performed?

Yes? or No?

Is the %RSD of the CFs for the two surrogates (tetrachloro-m-xylene and decachlorobiphenyl) less than or equal to 30.0%.

Yes? or No?

Action

- a. If the %RSD criteria are not met, qualify detects as estimated (J) and use professional judgment to qualify non-detected target compounds.
- b. If the %RSD criteria are within allowable limits, no qualification of the data is necessary

Continuing Calibration Checks

Criteria

Is the continuing calibration standard analyzed at the acceptable time intervals? Yes? or No?

Action

- a. If more than 14 hours has elapsed from the injection of the instrument blank that begins an analytical sequence (opening CCV) and the injection of either a PEM or mid-point concentration of the Individual Standard Mixtures (A and B) or (C), qualify all data as unusable (R).
- b. If more than 12 hours has elapsed from the injection of the instrument blank that begins an analytical sequence (opening CCV) and the injection of the last sample or blank that is part of the same analytical sequence, qualify all data as unusable (R).
- c. If more than 72 hours has elapsed from the injection of the sample with a Toxaphene detection and the Toxaphene Calibration Verification Standard (CS3), qualify all data as unusable (R).

Criteria

Is the Percent Difference (%D) within ±25.0% for the PEM sample?

Yes? or No?

Action

a. Qualify associated detects as estimated (J) and non-detects as estimated (UJ).

Criteria

For the Calibration Verification Standard (CS3); is the Percent Difference (%D) within ± 25.0%? Yes? or No?

Action

Qualify associated detects as estimated (UJ) and non-detects as estimated (UJ).

Criteria

Is the PEM 4,4'-DDT % Breakdown >20.0% and 4,4'-DDT is detected?

Yes? or No?

Action

- a. Qualify detects for 4,4'-DDT; detects for 4,4'-DDD; and detects for 4,4'-DDE as estimated (J)
- b. Non-detected associated compounds are not qualified

Criteria

Is the PEM 4,4'-DDT % Breakdown >20.0% and 4,4'-DDT is not detected

Yes? or No?

Action

- a. Qualify non-detects for 4,4'- DDT as unusable (R)
- b. Qualify detects for 4,4'-DDD as tentatively identified (NJ)
- c. Qualify detects for 4,4'-DDE as tentatively identified (NJ)

Criteria

Is the PEM Endrin % Breakdown >20.0% and Endrin is detected?

Yes? or No?

Action

- a. Qualify detects for Endrin; detects for Endrin aldehyde; and detects for Endrin ketone as estimated (J)
- b. Non-detected associated compounds are not qualified

Criteria

Is the PEM Endrin % Breakdown >20.0% and Endrin is not detected

Yes? or No?

Action

- a. Qualify non-detects for Endrin as unusable (R)
- b. Qualify detects for Endrin aldehyde as tentatively identified (NJ)
- c. Qualify detects for Endrin ketone as tentatively identified (NJ)

A separate worksheet should be filled for each initial curve

All criteria were metX
Criteria were not met
and/or see below

BLANK ANALYSIS RESULTS (Sections 1 & 2)

The assessment of the blank analysis results is to determine the existence and magnitude of contamination problems. The criteria for evaluation of blanks apply only to blanks associated with the samples, including trip, equipment, and laboratory blanks. If problems with any blanks exist, all data associated with the case must be carefully evaluated to determine whether or not there is an inherent variability in the data for the case, or if the problem is an isolated occurrence not affecting other data.

List the contam	ination in the bla	anks below. Hig	h and low levels blanks	must be treated separately.
CRQL concentr	rationN	/A		
Laboratory blan	ıks			
DATE Analyzed	LABID	LEVEL/ MATRIX	COMPOUND	CONCENTRATION UNITS
				nit_of_0.01_and_0.25_ug/L_
Field/Equipmen	t/Trip blank			
DATE Analyzed	LAB ID	LEVEL/ MATRIX	COMPOUND	CONCENTRATION UNITS
_No_field/trip/ed	quipment_blank	s_analyzed_wit	h_this_data_package	

All criteria were metX
Criteria were not met
and/or see below

BLANK ANALYSIS RESULTS (Section 3)

Blank Actions

Action Levels (ALs) should be based upon the highest concentration of contaminant determined in any blank. Do not qualify any blank with another blank. The ALs for samples which have been diluted should be corrected for the sample dilution factor and/or % moisture, where applicable. No positive sample results should be reported unless the concentration of the compound in the samples exceeds the ALs:

The concentration of non-target compounds in all blanks must be less than or equal to 10 μ g/L. The concentration of each target compound found in the method or field blanks must be less than its CRQL listed in the method.

Data concerning the field blanks are not evaluated as part of the CCS process. If field blanks are present, the data reviewer should evaluate this data in a similar fashion as the method blanks.

Specific actions are as follows:

Blank Actions for Pesticide Analyses

Blank Type	Blank Result	Sample Result	Action for Samples
	Detects	Not detected	No qualification required
Method, Sulfur Cleanup,	< CRQL	< CRQL	Report CRQL value with a U
		≥CRQL	No qualification required
		< CRQL	Report CRQL value with a U
		≥ CRQL and ≤ blank	Report blank value for
Instrument, Field, TCLP/SPLP	> CRQL	concentration	sample concentration with a U
		≥ CRQL and > blank concentration	No qualification required
	= CRQL	≤CRQL	Report CRQL value with a U
		> CRQL	No qualification required
	Gross contamination	Detects	Report blank value for sample concentration with a U

All criteria were metX
Criteria were not met
and/or see below

CONTAMINATION SOURCE/LEVEL	COMPOUND	CONC/UNITS	AL/UNITS	SQL	AFFECTED SAMPLES
		 			

All criteria were met __X__ Criteria were not met and/or see below ____

SURROGATE SPIKE RECOVERIES

Laboratory performance of individual samples is established by evaluation of surrogate spike recoveries. All samples are spiked with surrogate compounds prior to sample analysis. The accuracy of the analysis is measured by the surrogate percent recovery. Since the effects of the sample matrix are frequently outside the control of the laboratory and may present relatively unique problems, the validation of data is frequently subjective and demands analytical experience and professional judgment.

List the percent recoveries (%Rs) which do not meet the criteria for surrogate recovery.

Matrix:_Aqueou	IS				
Lab	Lab				
Sample ID	File ID	S1 a	S1 b	S2 a	S2 b
JC22613-7	6G36640.D	76	77	83	91
JC22613-8	6G36643.D	91	90	88	96
OP95079-BS1	6G36636.D	76	75	45	48
OP95079-MB1	6G36635.D	83	82	55	58
OP95079-MS	6G36641.D	98	98	85	94
OP95079-MSD	6G36642.D	73	72	68	74
Surrogate Com	pounds		Recove	ery Limit	S
S1 = Tetrachloro-m-xylene		26-132	%		
S2 = Decachlor	obiphenyl		10-118	1%	
(a) Recovery fro	om GC signal #1				(b) Recovery from GC signal i

Note: Surrogate recoveries within laboratory control limits.

Actions:

- a. For any surrogate recovery greater than 150%, qualify detected target compounds as biased high (J+).
- b. Do not qualify non-detected target compounds for surrogate recovery > 150 %.
- c. If both surrogate recoveries are greater than or equal to 30% and less than or equal to 150%, no qualification of the data is necessary.
- d. For any surrogate recovery greater than or equal to 10% and less than 30%, qualify detected target compounds as biased low (J-).
- e. For any surrogate recovery greater than or equal to 10% and less than 30%, qualify non-detected target compounds as approximated (UJ).
- f. If low surrogate recoveries are from sample dilution, professional judgment should be used to determine if the resulting data should be qualified. If sample dilution is not a factor:
 - i. Qualify detected target compounds as biased low (J-).
 - ii. Qualify non-detected target compounds as unusable (R).
- g. If surrogate RTs in PEMs, Individual Standard Mixtures, samples, and blanks are outside of the RT Windows, the reviewer must use professional judgment to qualify data.

- h. If surrogate RTs are within RT windows, no qualification of the data is necessary.
- i. If the two surrogates were not added to all samples, MS/MSDs, standards, LCSs, and blanks, use professional judgment in qualifying data as missing surrogate analyte may not directly apply to target analytes.

Summary Surrogate Actions for Pesticide Analyses

	Action*			
Criteria	Detected Target	Non-detected Target		
	Compounds	Compounds		
%R > 150%	J+	No qualification		
30% < %R < 150%	No qualification			
10% < %R < 30%	J-	UJ		
%R < 10% (sample dilution not a factor)	J-	R		
%R < 10% (sample dilution is a factor)	Use professional judgment			
RT out of RT window	Use professional judgment			
RT within RT window	No qualification			

^{*} Use professional judgment in qualifying data, as surrogate recovery problems may not directly apply to target analytes.

Il criteria were metX
Criteria were not met
and/or see below

MATRIX SPIKE/MATRIX SPIKE DUPLICATE (MS/MSD)

This data is generated to determine long term precision and accuracy in the analytical method for various matrices. This data alone cannot be used to evaluate the precision and accuracy of individual samples. If any % R in the MS or MSD falls outside the designated range, the reviewer should determine if there are matrix effects, i.e. LCS data are within the QC limits but MS/MSD data are outside QC limit.

MS/MSD Recoveries and Precision Criteria

laboratory control limits.

Data for MS and MSDs will not be present unless requested by the Region.

Notify the Contract Laboratory Program Project Officer (CLP PO) if a field blank was used for the MS and MSD, unless designated as such by the Region.

NOTE: For a Matrix Spike that does not meet criteria, apply the action to only the field sample used to prepare the Matrix Spike sample. If it is clearly stated in the data validation materials that the samples were taken through incremental sampling or some other method guaranteeing the homogeneity of the sample group, then the entire sample group may be qualified.

Sample ID:	JC22613-7MS/MSD_			Matrix/	Level:Groundwater
MS OR MSD	COMPOUND	% R	RPD	QC LIMITS	ACTION

Action

No qualification of the data is necessary on MS and MSD data alone. However, using professional judgment, the validator may use the MS and MSD results in conjunction with other QC criteria and determine the need for some qualification of the data.

Note: MS/MSD sample analyzed with this data package. % recoveries and RPD within

A separate worksheet should be used for each MS/MSD pair.

All criteria were metX
Criteria were not met
and/or see below

LABORATORY CONTROL SAMPLE (LCS) ANALYSIS

This data is generated to determine accuracy of the analytical method for various matrices.

1. LCS Recoveries Criteria

LCS Spike Compound	Recovery Limits (%)
gamma-BHC	50 – 120
Heptachlor epoxide	50 – 150
Dieldrin	30 – 130
4,4'-DDE	50 – 150
Endrin	50 – 120
Endosulfan sulfate	50 – 120
trans-Chlordane	30 – 130
Tetrachloro-m-xylene (surrogate)	30 – 150
Decachlorobiphenyl (surrogate)	30 – 150

	:0.25_ug/l; which do not meet the criteria		· · · · · · · · · · · · · · · · · · ·
LCS ID	COMPOUND	% R	QC LIMIT

Action

The following guidance is suggested for qualifying sample data for which the associated LCS does not meet the required criteria.

- a. If the LCS recovery exceeds the upper acceptance limit, qualify detected target compounds as estimated (J). Do not qualify non-detected target compounds.
- b. If the LCS recovery is less than the lower acceptance limit, qualify detected target compounds as estimated (J) and non-detects as unusable (R).
- c. Use professional judgment to qualify data for compounds other than those compounds that are included in the LCS.
- d. Use professional judgment to qualify non-LCS compounds. Take into account the compound class, compound recovery efficiency, analytical problems associated with each compound, and comparability in the performance of the LCS compound to the non-LCS compound.
- e. If the LCS recovery is within allowable limits, no qualification of the data is necessary.

2. Frequency Criteria:

Where LCS analyzed at the required frequency and for each matrix? <u>Yes</u> or No. If no, the data may be affected. Use professional judgment to determine the severity of the effect and qualify data accordingly. Discuss any actions below and list the samples affected.

Note: Blank spike analyzed for aqueous matrix. % recoveries within laboratory control limits. Recovery for gamma-chlordane obtained from second column, first column used for confirmation only.

All criteria were met
Criteria were not met
and/or see belowN/A

FLORISIL CARTRIDGE PERFORMANCE CHECK

NOTE: Florisil cartridge cleanup is mandatory for all extracts.

Criteria

Is the Florisil cartridge performance check conducted at least once on each lot of cartridges used for sample cleanup or every 6 months, whichever is most frequent?

Yes? or No?

N/A

Criteria

Are the results for the Florisil Cartridge Performance Check solution included with the data package?

Yes? or No? N/A

Note: If % criteria are not met, examine the raw data for the presence of polar interferences and use professional judgment in qualifying the data as follows:

Action:

- a. If the Percent Recovery is greater than 120% for any of the pesticide target compounds in the Florisil Cartridge Performance Check, qualify detected compounds as estimated (J). Do not qualify non-detected target compounds.
- b. If the Percent Recovery is greater than or equal to 80% and less than or equal to 120% for all the pesticide target compounds, no qualification of the data is necessary.
- c. If the Percent Recovery is greater than or equal to 10% and less than 80% for any of the pesticide target compounds in the Florisil Cartridge Performance Check, qualify detected target compounds as estimated (J) and non-detected target compounds as approximated (UJ).
- d. If the Percent Recovery is less than 10% for any of the pesticide target compounds in the Florisil Cartridge Performance Check, qualify detected compounds as estimated (J) and qualify non-detected target compounds as unusable (R).
- e. If the Percent Recovery of 2,4,5-trichlorophenol in the Florisil Cartridge Performance Check is greater than or equal to 5%, use professional judgment to qualify detected and non-detected target compounds, considering interference on the sample chromatogram.

Note: State in the Data Review Narrative potential effects on the sample data resulting from the Florisil Cartridge Performance Check analysis not yielding acceptable results.

Note: No information for florisil cartridge performance check included in data package. There is evidence tahtFlorisil cartridge was used for sample extraction/clean-up. No qualification of the data performed, professional judgment.

All criteria were metNA	<u> </u>
Criteria were not met	
and/or see below	

GEL PERMEATION CHROMATOGRAPHY (GPC) PERFORMANCE CHECK

NOTE: GPC cleanup is mandatory for all soil samples.

If GPC criteria are not met, examine the raw data for the presence of high molecular weight contaminants; examine subsequent sample data for unusual peaks; and use professional judgment in qualifying the data. Notify the Contract Laboratory Program Project Officer (CLP PO) if the laboratory chooses to analyze samples under unacceptable GPC criteria.

Action:

- a. If the Percent Recovery is less than 10% for the pesticide compounds and surrogates during the GPC calibration check, the non-detected target compounds may be suspect, qualify detected compounds as estimated (J).
- b. If the Percent Recovery is less than 10% for the pesticide compounds and surrogates during the GPC calibration check, qualify all non-detected target compounds as unusable (R).
- c. If the Percent Recovery is greater than or equal to 10% and is less than 80% for any of the pesticide target compounds in the GPC calibration, qualify detected target compounds as estimated (J) and non-detected target compounds as approximated (UJ).
- d. If the Percent Recovery is greater than or equal to 80% and less than or equal to 120% for all the pesticide target compounds, no qualification of the data is necessary.
- e. If high recoveries (i.e., greater than 120%) were obtained for the pesticides and surrogates during the GPC calibration check, qualify detected compounds as estimated (J). Do not qualify non-detected target compounds.

Note: State in the Data Review Narrative potential effects on the sample data resulting from the GPC cleanup analyses not yielding acceptable results.

Note: No information for performance of GPC cleanup included in data package. No qualification of the data performed, professional judgment.

All criteria were metX	_
Criteria were not met	
and/or see below	_

TARGET COMPOUND IDENTIFICATION

Criteria:

- 1. Is Retention Times (RTs) of both of the surrogates and reported target compounds in each sample within the calculated RT Windows on both columns?

 Yes? or No?
- 2. Is the Tetrachloro-m-xylene (TCX) RT ±0.05 minutes of the Mean RT (RT) determined from the initial calibration and Decachlorobiphenyl (DCB) within ±0.10 minutes of the RT determined from the initial calibration?

 Yes? or No?
- 3. Is the Percent Difference (%D) for the detected mean concentrations of a pesticide target compound between the two Gas Chromatograph (GC) columns within the inclusive range of ± 25.0 %?

 Yes? or No?
- 4. When no analytes are identified in a sample; are the chromatograms from the analyses of the sample extract and the low-point standard of the initial calibration associated with those analyses on the same scaling factor?

 Yes? or No?
- 5. Does the chromatograms display the Single Component Pesticides (SCPs) detected in the sample and the largest peak of any multi-component analyte detected in the sample at less than full scale.

 Yes? or No?
- 6. If an extract is diluted; does the chromatogram display SCPs peaks between 10-100% of full scale, and multi-component analytes between 25-100% of full scale? Yes? or No? N/A
- 7. For any sample; does the baseline of the chromatogram return to below 50% of full scale before the elution time of alpha-BHC, and also return to below 25% of full scale after the elution time of alpha-BHC and before the elution time of DCB?

 Yes? or No?
- 8. If a chromatogram is replotted electronically to meet these requirements; is the scaling factor used displayed on the chromatogram, and both the initial chromatogram and the replotted chromatogram submitted in the data package.

 Yes? or No?

Action:

- a. If the qualitative criteria for both columns were not met, all target compounds that are reported as detected should be considered non-detected.
- b. Use professional judgment to assign an appropriate quantitation limit using the following guidance:
 - If the detected target compound peak was sufficiently outside the pesticide RT Window, the reported values may be a false positive and should be replaced with the sample Contract Required Quantitation Limits (CRQL) value.

- ii. If the detected target compound peak poses an interference with potential detection of another target peak, the reported value should be considered and qualified as unusable (R).
- c. If the data reviewer identifies a peak in both GC column analyses that falls within the appropriate RT Windows, but was reported as a non-detect, the compound may be a false negative. Use professional judgment to decide if the compound should be included.

Note: State in the Data Review Narrative all conclusions made regarding target compound identification.

- d. If the Toxaphene peak RT windows determined from the calibration overlap with SCPs or chromatographic interferences, use professional judgment to qualify the data.
- e. If target compounds were detected on both GC columns, and the Percent Difference between the two results is greater than 25.0%, consider the potential for coelution and use professional judgment to decide whether a much larger concentration obtained on one column versus the other indicates the presence of an interfering compound. If an interfering compound is indicated, use professional judgment to determine how best to report, and if necessary, qualify the data according to these guidelines.
- f. If Toxaphene exhibits a marginal pattern-matching quality, use professional judgment to establish whether the differences are due to environmental "weathering" (i.e., degradation of the earlier eluting peaks relative to the later eluting peaks). If the presence of Toxaphene is strongly suggested, report results as presumptively present (N).

GAS CHROMATOGRAPH/MASS SPECTROMETER (GC/MS) CONFIRMATION

NOTE: This confirmation is not usually provided by the laboratory. In cases where it is provided, use professional judgment to determine if data qualified with "C" can be salvaged if it was previously qualified as unusable (R).

Action:

- a. If the quantitative criteria for both columns were met (≥ 5.0 ng/µL for SCPs and ≥ 125 ng/µL for Toxaphene), determine whether GC/MS confirmation was performed. If it was performed, qualify the data using the following guidance:
 - i. If GC/MS confirmation was not required because the quantitative criteria for both columns was not met, but it was still performed, use professional judgment when evaluating the data to decide whether the detect should be qualified with "C".
 - ii. If GC/MS confirmation was performed, but unsuccessful for a target compound detected by GC/ECD analysis, qualify those detects as "X".

All criteria were metX
Criteria were not met
and/or see below

COMPOUND QUANTITATION AND REPORTED CONTRACT REQUIRED QUANTITATION LIMITS (CRQLS)

The sample quantitation evaluation is to verify laboratory quantitation results. In the space below, please show a minimum of one sample calculation:

JC22613-7MS

alpha-BHC

RF = 1.248

[] =

(76547368)(50)/(121.2 X 10⁶)(1.248)

= 25.30 ppb

Ok

Action:

- a. If sample quantitation is different from the reported value, qualify result as unusable (R).
- b. When a sample is analyzed at more than one dilution, the lowest CRQLs are used unless a QC exceedance dictates the use of the higher CRQLs from the diluted sample.
- c. Replace concentrations that exceed the calibration range in the original analysis by crossing out the "E" and its corresponding value on the original reporting form and substituting the data from the diluted sample.
- d. Results between the MDL and CRQL should be qualified as estimated (J).
- e. Results less than the MDL should be reported at the CRQL and qualified (U). MDLs themselves are not reported.
- f. For non-aqueous samples, if the percent moisture is less than 70.0%, no qualification of the data is necessary. If the percent moisture is greater than or equal to 70.0% and less than 90.0%, qualify detects as estimated (J) and non-detects as approximated (UJ). If the percent moisture is greater than or equal to 90.0%, qualify detects as estimated (J) and non-detects as unusable (R) (see Table).

Percent Moisture Actions for Pesticide Analysis for Non-Aqueous Samples

Criteria	Action		
	Detected Associated Compounds	Non-detected Associated Compounds	
% Moisture < 70.0	No qualification		
70.0 < % Moisture < 90.0	J	้ เก	
% Moisture > 90.0	J	R	

samples which have $\leq 50 \%$	0 301103		

Note: If any discrepancies are found, the Region's designated representative may contact the laboratory to obtain additional information that could resolve any differences. If a discrepancy remains unresolved, the reviewer must use professional judgment to decide which value is the most accurate. Under these circumstances, the reviewer may determine that qualification of data is warranted. Note in the Data Review Narrative a description of the reasons for data qualification and the qualification that is applied to the data.

Dilution performed

SAMPLE ID	DILUTION FACTOR	REASON FOR DILUTION
	-	
		+

All criteria were met _	_NVA
Criteria were not met	
and/or see below	

FIELD DUPLICATE PRECISION

NOTE: In the absence of QAPP guidance for validating data from field duplicates, the following action will be taken.

Field duplicates samples may be taken and analyzed as an indication of overall precision. These analyses measure both field and lab precision; therefore, the results may have more variability than laboratory duplicates which only laboratory performance. It is also expected that soil duplicate results will have a greater variance than water matrices due to difficulties associated with collecting identical field duplicate samples. Identify which samples within the data package are field duplicates. Estimate the relative percent difference (RPD) between the values for each compound. If large RPDs (> 50%) is observed, confirm identification of samples and note difference in the executive summary.

Sample IDs				/latrix:	•
COMPOUND	SQL ug/L	SAMPLE CONC.	DUPLICATE CONC.	RPD	ACTION
			data package. MS/I in the required criter		coveries RPD used to %.

Actions:

- a. Qualify as estimated positive results (J) and nondetects (UJ) for the compound that exceeded the above criteria. For organics, only the sample and duplicate will be qualified.
- b. If an RPD cannot be calculated because one or both of the sample results is not detected, the following actions apply:
 - i. If one sample result is not detected and the other is greater than 5x the SQL qualify (J/UJ).
 - ii. If one sample value is not detected and the other is greater than 5x the SQL and the SQLs for the sample and duplicate are significantly different, use professional judgment to determine if qualification is appropriate.
 - iii. If one sample value is not detected and the other is less than 5x, use professional judgment to determine if qualification is appropriate.
 - iv. If both sample and duplicate results are not detected, no action is needed.

OVERALL ASSESSMENT OF DATA

Action:

- 1. Use professional judgment to determine if there is any need to qualify data which were not qualified based on the Quality Control (QC) criteria previously discussed.
- 2. Write a brief narrative to give the user an indication of the analytical limitations of the data.

Note: The Contract Laboratory Program Project Officer (CLP PO) must be informed if any inconsistency of the data with the Sample Delivery Group (SDG) Narrative. If sufficient information on the intended use and required quality of the data is available, the reviewer should include their assessment of the usability of the data within the given context. This may be used as part of a formal Data Quality Assessment (DQA).

Overall assessment of the data:

Results are valid; the data can be used for decision making purposes.