

MEMORANDUM

checked by #
09/28/15

TO: Mr. Terry Taylor
Anderson, Mulholland and Associates

DATE: September 16, 2015

FROM: R. Infante

FILE: JC2870

RE: Data Validation
Air samples
BMSMC, Building 5 Area, PR
SM04.00.06
Building 8 Vapor Intrusions

SUMMARY

Full validation was performed on the data for several gas samples analyzed for selected volatile organic compounds by method Compendium Method TO-15: Determination Of Volatile Organic Compounds (VOCs) In Air Collected In Specially-Prepared Canisters And Analyzed By Gas Chromatography/Mass Spectrometry (GC/MS), January, 1999 and by method TO-3: METHOD FOR THE DETERMINATION OF VOLATILE ORGANIC COMPOUNDS IN AMBIENT AIR USING CRYOGENIC PRECONCENTRATION TECHNIQUES AND GAS CHROMATOGRAPHY WITH FLAME IONIZATION AND ELECTRON CAPTURE DETECTION, April, 1984. The samples were collected at the Bristol Myer Squib-Building 5 area, Humacao, PR site on August 31, 2015 and submitted to Accutest Laboratories, of Dayton, New Jersey that analyzed and reported the results under delivery groups (SDG) JC2870.

The sample results were assessed according to USEPA data validation guidance documents in the following order of precedence: Compendium Method TO-15. Determination Of Volatile Organic Compounds (VOCs) In Air Collected In Specially-Prepared Canisters And Analyzed By Gas Chromatography/Mass Spectrometry (GC/MS), January, 1999; Validating Air Samples. Volatile Organic Analysis of Ambient Air in Canisters by Method TO-15, (SOP # HW-31. Revision #4. October, 2006. The QC criteria and data validation actions listed on the data review worksheets are from the primary guidance document, unless otherwise noted. In general the data is valid as reported and may be used for decision making purposes.

The data results are acceptable for use.

SAMPLES

The samples included in the review are listed below

| Client Sample ID | Lab. Sample ID | Collected Date | Matrix | Analysis |
|---------------------|----------------|-------------------|--------|-------------|
| B8-SSV2/08312015 | JC2870-1 | 08/31/2015 | Air | TO-3 |
| B8-SSV2/08312015 | JC2870-1A | 08/31/2015 | Air | TO-15 |
| B8-SSV1/08312015 | JC2870-2 | 08/31/2015 | Air | TO-15; TO-3 |
| B8-SSV1D/08312015 | JC2870-3 | 08/31/2015 | Air | TO-15; TO-3 |

REVIEW ELEMENTS

Sample data were reviewed for the following parameters, where applicable to the method

- Agreement of analysis conducted with chain of custody (COC) form
- Holding time and sample preservation
- Gas chromatography/mass spectrometry (GC/MS) tunes
- Initial and continuing calibrations
- Method blanks/trip blanks/field blank
- Canister cleaning certification criteria
- Surrogate spike recovery
- Internal standard performance and retention times
- Field duplicate results
- Laboratory control sample/laboratory control sample duplicate (LCS/LCSD) results
- Quantitation limits and sample results

DISCUSSION

Agreement of Analysis Conducted with COC Request

Sample reports corresponded to the analytical request designated on the chain-of-custody form.

Holding Times and Sample Preservation

Sample preservation was acceptable.

Samples analyzed within method recommended holding time.

GC/MS Tunes

The frequency and abundance of bromofluorobenzene (BFB) tunes were within the QC acceptance criteria. All samples were analyzed within the tuning criteria associated with the method.

Initial and Continuing Calibrations

VOCs (Method TO-15)

The percent relative standard deviations (%RSDs) and response factors (RFs) of all target analytes were within the QC acceptance criteria in the initial calibration. Correlation coefficients (r^2) of target analytes were within the QC acceptance criteria. Ongoing accuracy of the instrument was determined by the analysis of a continuing calibration standard. All target analytes meet the method calibration/continued.

Method Blank/Trip Blank/Field Blank

Target analytes were not detected in laboratory method blanks.

Summa canister met cleaning certification criteria.

No trip/field blank analyzed with this data package.

Surrogate Spike Recovery

The surrogate recoveries were within the laboratory QC acceptance limits in all samples analyzed except for the followings:

TO-3:

- JC2870-1:- 4-BFB outside laboratory control limits (137 % recovery); control limits 61 – 137 %.
- JC2870-1D:- 4-BFB outside laboratory control limits (137 % recovery); control limits 61 – 137 %.

No action taken; % recoveries within control limits in second column. % recoveries within generally acceptable control limits.

TO-15:

- JC2870-2:- 4-BFB outside laboratory control limits (39 % recovery); control limits 65 – 128 %.
- JC2870-2:- 4-BFB outside laboratory control limits (36 % recovery); control limits 65 – 128 %.
- JC2870-3:- 4-BFB outside laboratory control limits (55 % recovery); control limits 65 – 128 %.
- JC2870-1A:- 4-BFB outside laboratory control limits (22 % recovery); control limits 65 – 128 %.
- JC2870-1A:- 4-BFB outside laboratory control limits (23 % recovery); control limits 65 – 128 %.

No action taken, professional judgment.

Internal Standard Performance

VOCs (TO-15)

Samples were spiked with the method specified internal standard. Internal standard are performance and retention times met the QC acceptance criteria in all sample analyses and calibration standards except for the followings:

| DATE | SAMPLE ID | IS OUT | IS AREA | ACCEPTABLE RANGE | ACTION |
|----------|-----------|--------------------|---------|------------------|-----------|
| 09/03/15 | JC2870-1A | 5-Chlorobenzene-d5 | 870406 | 150209 - 350489 | No action |
| 09/03/15 | JC2870-2 | 5-Chlorobenzene-d5 | 588889 | 150209 - 350489 | No action |

No action taken, professional judgment.

Laboratory/Field Duplicate Results

TO-3

Laboratory duplicates (JC2870-1/-1 DUP and JC3132-1/-1 DUP) were analyzed as part of this data set. Target analytes meet the RPD performance criteria and generally acceptable laboratory control limits for analytes 5 x SQL. Field duplicates were samples B8-SSV108312015/B8-SSV1D/08312015,

target analytes meet the RPD performance criteria and generally acceptable control limits for analytes 5 x SQL.

TO-15

Laboratory duplicates (JC2480-6/-6 Duplicate and JC1343-1/-1 Duplicate) were analyzed as part of this data set. Target analytes meet the RPD performance criteria and generally acceptable laboratory control limits for analytes 5 x SQL. Field duplicates were samples B8-SSV108312015/B8-SSV1D/08312015, target analytes meet the RPD performance criteria and generally acceptable control limits for analytes 5 x SQL.

LCS/LCSD Results

VOCs

LCS/LCSD (blank spike) associated with this data package were analyzed by the laboratory. Recoveries and RPD within laboratory control limits.

Quantitation Limits and Sample Results

TO-3:

Dilution was performed on sample JC2870-1 due to the presence of high level target species. All other samples diluted by a factor of less than 3.

Detected results for Methane in sample JC2870-1 exceed the instrument calibration range and are considered estimated values.


To-15

All samples diluted by a factor of less than 3.

Calculations were spot checked.

Certification

The following samples JC2870-1; JC2870-1A; JC2870-2; and JC2870-3 were analyzed following standard procedures accepted by regulatory agencies. The quality control requirements met the methods criteria except in the occasions described in this document. The results are valid. Some of the results were qualified.


Rafael Infante
Chemist License 1888



Accutest LabLink@872636 09:45 14-Sep-2015

Report of Analysis

Page 1 of 1

| | | | |
|-------------------|--|-----------------|----------|
| Client Sample ID: | B8-SSV2/08312015 | Date Sampled: | 08/31/15 |
| Lab Sample ID: | JC2870-1 | Date Received: | 09/02/15 |
| Matrix: | AIR - Soil Vapor Comp. Summa ID: A1216 | Percent Solids: | n/a |
| Method: | EPA TO-3 | | |
| Project: | BMSMC, Building 5 Area, PR | | |

| Run # | File ID | DF | Analyzed | By | Prep Date | Prep Batch | Analytical Batch |
|--------|-----------|------|----------|-----|-----------|------------|------------------|
| Run #1 | QT19182.D | 2.96 | 09/03/15 | TCH | n/a | n/a | GQT1156 |
| Run #2 | | | | | | | |

| Run # | Initial Volume |
|--------|----------------|
| Run #1 | 0.50 ml |
| Run #2 | |

| CAS No. | MW | Compound | Result | RL | MDL | Units | Q | Result | RL | MDL | Units |
|----------|----|----------|--------|-----|------|-------|---|--------|-----|------|-------|
| 74-82-8 | 16 | Methane | 385000 | 15 | 0.99 | ppmv | E | 252000 | 9.8 | 0.65 | mg/m3 |
| 74-98-6 | 44 | Propane | ND | 1.5 | 0.12 | ppmv | | ND | 2.7 | 0.22 | mg/m3 |
| 106-97-8 | 58 | n-Butane | ND | 1.5 | 0.18 | ppmv | | ND | 3.6 | 0.43 | mg/m3 |

| CAS No. | Surrogate Recoveries | Run# 1 | Run# 2 | Limits |
|----------|----------------------|-------------------|--------|---------|
| 460-00-4 | 4-Bromofluorobenzene | 137% ^a | | 61-137% |
| 460-00-4 | 4-Bromofluorobenzene | 111% | | 61-137% |

(a) Outside control limits due to matrix interference, confirmed by sample duplicate.



ND = Not detected MDL = Method Detection Limit
 RL = Reporting Limit
 E = Indicates value exceeds calibration range

J = Indicates an estimated value
 B = Indicates analyte found in associated method blank
 N = Indicates presumptive evidence of a compound

Accutest LabLink@872636 09:45 14-Sep-2015

Report of Analysis

Page 1 of 3

| | | | |
|-------------------|---|-----------------|----------|
| Client Sample ID: | B8-SSV2/08312015 | Date Sampled: | 08/31/15 |
| Lab Sample ID: | JC2870-1A | Date Received: | 09/02/15 |
| Matrix: | AIR - Soil Vapor Comp. Summa ID: A1216,A681 | Percent Solids: | n/a |
| Method: | TO-15 | | |
| Project: | BMSMC, Building 5 Area, PR | | |

| Run # | File ID | DF | Analyzed | By | Prep Date | Prep Batch | Analytical Batch |
|---------------------|-----------|------|----------|-----|-----------|------------|------------------|
| Run #1 | 3W49835.D | 59.2 | 09/02/15 | YMH | n/a | n/a | V3W1889 |
| Run #2 ^a | 3W49845.D | 59.2 | 09/03/15 | YMH | n/a | n/a | V3W1890 |


| | Initial Volume |
|--------|----------------|
| Run #1 | 200 ml |
| Run #2 | 200 ml |

VOA TO15 List

| CAS No. | MW | Compound | Result | RL | MDL | Units | Q | Result | RL | MDL | Units |
|------------|-------|----------------------------|--------|----|-----|-------|------|--------|-----|-------|-------|
| 67-64-1 | 58.08 | Acetone | 91.3 | 24 | 3.8 | ppbv | 217 | 57 | 9.0 | ug/m3 | |
| 106-99-0 | 54.09 | 1,3-Butadiene | ND | 24 | 3.7 | ppbv | ND | 53 | 8.2 | ug/m3 | |
| 71-43-2 | 78.11 | Benzene | ND | 24 | 3.5 | ppbv | ND | 77 | 11 | ug/m3 | |
| 75-27-4 | 163.8 | Bromodichloromethane | ND | 24 | 3.8 | ppbv | ND | 160 | 25 | ug/m3 | |
| 75-25-2 | 252.8 | Bromoform | ND | 24 | 2.4 | ppbv | ND | 250 | 25 | ug/m3 | |
| 74-83-9 | 94.94 | Bromomethane | ND | 24 | 2.6 | ppbv | ND | 93 | 10 | ug/m3 | |
| 593-60-2 | 106.9 | Bromoethene | ND | 24 | 2.4 | ppbv | ND | 100 | 10 | ug/m3 | |
| 100-44-7 | 126 | Benzyl Chloride | ND | 24 | 3.1 | ppbv | ND | 120 | 16 | ug/m3 | |
| 75-15-0 | 76.14 | Carbon disulfide | ND | 24 | 3.4 | ppbv | ND | 75 | 11 | ug/m3 | |
| 108-90-7 | 112.6 | Chlorobenzene | ND | 24 | 3.8 | ppbv | ND | 110 | 18 | ug/m3 | |
| 75-00-3 | 64.52 | Chloroethane | ND | 24 | 2.5 | ppbv | ND | 63 | 6.6 | ug/m3 | |
| 67-66-3 | 119.4 | Chloroform | ND | 24 | 3.7 | ppbv | ND | 120 | 18 | ug/m3 | |
| 74-87-3 | 50.49 | Chloromethane | ND | 24 | 3.4 | ppbv | ND | 50 | 7.0 | ug/m3 | |
| 107-05-1 | 76.53 | 3-Chloropropene | ND | 24 | 3.3 | ppbv | ND | 75 | 10 | ug/m3 | |
| 95-49-8 | 126.6 | 2-Chlorotoluene | ND | 24 | 3.9 | ppbv | ND | 120 | 20 | ug/m3 | |
| 56-23-5 | 153.8 | Carbon tetrachloride | ND | 24 | 2.9 | ppbv | ND | 150 | 18 | ug/m3 | |
| 110-82-7 | 84.16 | Cyclohexane | 1530 | 24 | 3.8 | ppbv | 5270 | 83 | 13 | ug/m3 | |
| 75-34-3 | 98.96 | 1,1-Dichloroethane | ND | 24 | 3.7 | ppbv | ND | 97 | 15 | ug/m3 | |
| 75-35-4 | 96.94 | 1,1-Dichloroethylene | ND | 24 | 3.4 | ppbv | ND | 95 | 13 | ug/m3 | |
| 106-93-4 | 187.9 | 1,2-Dibromoethane | ND | 24 | 4.2 | ppbv | ND | 180 | 32 | ug/m3 | |
| 107-06-2 | 98.96 | 1,2-Dichloroethane | ND | 24 | 3.1 | ppbv | ND | 97 | 13 | ug/m3 | |
| 78-87-5 | 113 | 1,2-Dichloropropane | ND | 24 | 6.0 | ppbv | ND | 110 | 28 | ug/m3 | |
| 123-91-1 | 88.12 | 1,4-Dioxane | ND | 24 | 7.5 | ppbv | ND | 86 | 27 | ug/m3 | |
| 75-71-8 | 120.9 | Dichlorodifluoromethane | ND | 24 | 4.3 | ppbv | ND | 120 | 21 | ug/m3 | |
| 124-48-1 | 208.3 | Dibromochloromethane | ND | 24 | 4.9 | ppbv | ND | 200 | 42 | ug/m3 | |
| 156-60-5 | 96.94 | trans-1,2-Dichloroethylene | ND | 24 | 2.4 | ppbv | ND | 95 | 9.5 | ug/m3 | |
| 156-59-2 | 96.94 | cis-1,2-Dichloroethylene | ND | 24 | 2.9 | ppbv | ND | 95 | 11 | ug/m3 | |
| 10061-01-5 | 111 | cis-1,3-Dichloropropene | ND | 24 | 4.2 | ppbv | ND | 110 | 19 | ug/m3 | |
| 541-73-1 | 147 | m-Dichlorobenzene | ND | 24 | 3.3 | ppbv | ND | 140 | 20 | ug/m3 | |
| 95-50-1 | 147 | o-Dichlorobenzene | ND | 24 | 3.6 | ppbv | ND | 140 | 22 | ug/m3 | |
| 106-46-7 | 147 | p-Dichlorobenzene | ND | 24 | 2.3 | ppbv | ND | 140 | 14 | ug/m3 | |
| 10061-02-6 | 111 | trans-1,3-Dichloropropene | ND | 24 | 2.4 | ppbv | ND | 110 | 11 | ug/m3 | |

ND = Not detected MDL = Method Detection Limit
 RL = Reporting Limit
 E = Indicates value exceeds calibration range

J = Indicates an estimated value
 B = Indicates analyte found in associated method blank
 N = Indicates presumptive evidence of a compound



Report of Analysis

Page 2 of 3

| | |
|--|--------------------------------|
| Client Sample ID: B8-SSV2/08312015 | Date Sampled: 08/31/15 |
| Lab Sample ID: JC2870-1A | Date Received: 09/02/15 |
| Matrix: AIR - Soil Vapor Comp. Summa ID: A1216,A681 | Percent Solids: n/a |
| Method: TO-15 | |
| Project: BMSMC, Building 5 Area, PR | |

VOA TO15 List

| CAS No. | MW | Compound | Result | RL | MDL | Units | Q | Result | RL | MDL | Units |
|-----------|--------|---------------------------|--------|-----|-----|-------|---|--------|-----|-----|-------|
| 64-17-5 | 46.07 | Ethanol | ND | 59 | 20 | ppbv | | ND | 110 | 38 | ug/m3 |
| 100-41-4 | 106.2 | Ethylbenzene | ND | 24 | 5.7 | ppbv | | ND | 100 | 25 | ug/m3 |
| 141-78-6 | 88 | Ethyl Acetate | ND | 24 | 7.5 | ppbv | | ND | 86 | 27 | ug/m3 |
| 622-96-8 | 120.2 | 4-Ethyltoluene | ND | 24 | 2.6 | ppbv | | ND | 120 | 13 | ug/m3 |
| 76-13-1 | 187.4 | Freon 113 | ND | 24 | 3.2 | ppbv | | ND | 180 | 25 | ug/m3 |
| 76-14-2 | 170.9 | Freon 114 | ND | 24 | 3.0 | ppbv | | ND | 170 | 21 | ug/m3 |
| 142-82-5 | 100.2 | Heptane | 1920 | 24 | 3.5 | ppbv | | 7870 | 98 | 14 | ug/m3 |
| 87-68-3 | 260.8 | Hexachlorobutadiene | ND | 24 | 3.9 | ppbv | | ND | 260 | 42 | ug/m3 |
| 110-54-3 | 86.17 | Hexane | 817 | 24 | 3.3 | ppbv | | 2880 | 85 | 12 | ug/m3 |
| 591-78-6 | 100 | 2-Hexanone | ND | 24 | 5.2 | ppbv | | ND | 98 | 21 | ug/m3 |
| 67-63-0 | 60.1 | Isopropyl Alcohol | ND | 24 | 14 | ppbv | | ND | 59 | 34 | ug/m3 |
| 75-09-2 | 84.94 | Methylene chloride | ND | 24 | 16 | ppbv | | ND | 83 | 56 | ug/m3 |
| 78-93-3 | 72.11 | Methyl ethyl ketone | ND | 24 | 5.8 | ppbv | | ND | 71 | 17 | ug/m3 |
| 108-10-1 | 100.2 | Methyl Isobutyl Ketone | ND | 24 | 3.2 | ppbv | | ND | 98 | 13 | ug/m3 |
| 1634-04-4 | 88.15 | Methyl Tert Butyl Ether | 531 | 24 | 3.1 | ppbv | | 1910 | 87 | 11 | ug/m3 |
| 80-62-6 | 100.12 | Methylmethacrylate | ND | 24 | 3.6 | ppbv | | ND | 98 | 15 | ug/m3 |
| 115-07-1 | 42 | Propylene | 250 | 59 | 9.6 | ppbv | | 429 | 100 | 16 | ug/m3 |
| 100-42-5 | 104.1 | Styrene | ND | 24 | 3.0 | ppbv | | ND | 100 | 13 | ug/m3 |
| 71-55-6 | 133.4 | 1,1,1-Trichloroethane | ND | 24 | 3.7 | ppbv | | ND | 130 | 20 | ug/m3 |
| 79-34-5 | 167.9 | 1,1,2,2-Tetrachloroethane | ND | 24 | 3.6 | ppbv | | ND | 160 | 25 | ug/m3 |
| 79-00-5 | 133.4 | 1,1,2-Trichloroethane | ND | 24 | 4.2 | ppbv | | ND | 130 | 23 | ug/m3 |
| 120-82-1 | 181.5 | 1,2,4-Trichlorobenzene | ND | 24 | 5.2 | ppbv | | ND | 180 | 39 | ug/m3 |
| 95-63-6 | 120.2 | 1,2,4-Trimethylbenzene | ND | 24 | 2.7 | ppbv | | ND | 120 | 13 | ug/m3 |
| 108-67-8 | 120.2 | 1,3,5-Trimethylbenzene | ND | 24 | 3.6 | ppbv | | ND | 120 | 18 | ug/m3 |
| 540-84-1 | 114.2 | 2,2,4-Trimethylpentane | ND | 24 | 2.5 | ppbv | | ND | 110 | 12 | ug/m3 |
| 75-65-0 | 74.12 | Tertiary Butyl Alcohol | 21.1 | 24 | 5.9 | ppbv | J | 64.0 | 73 | 18 | ug/m3 |
| 127-18-4 | 165.8 | Tetrachloroethylene | ND | 4.7 | 2.8 | ppbv | | ND | 32 | 19 | ug/m3 |
| 109-99-9 | 72.11 | Tetrahydrofuran | ND | 24 | 5.1 | ppbv | | ND | 71 | 15 | ug/m3 |
| 108-88-3 | 92.14 | Toluene | 12.3 | 24 | 2.4 | ppbv | J | 46.4 | 90 | 9.0 | ug/m3 |
| 79-01-6 | 131.4 | Trichloroethylene | ND | 4.7 | 3.0 | ppbv | | ND | 25 | 16 | ug/m3 |
| 75-69-4 | 137.4 | Trichlorofluoromethane | ND | 24 | 2.4 | ppbv | | ND | 130 | 13 | ug/m3 |
| 75-01-4 | 62.5 | Vinyl chloride | 42.4 | 24 | 3.8 | ppbv | | 108 | 61 | 9.7 | ug/m3 |
| 108-05-4 | 86 | Vinyl Acetate | ND | 24 | 6.5 | ppbv | | ND | 84 | 23 | ug/m3 |
| | 106.2 | m,p-Xylene | 30.4 | 24 | 5.1 | ppbv | | 132 | 100 | 22 | ug/m3 |
| 95-47-6 | 106.2 | o-Xylene | ND | 24 | 3.0 | ppbv | | ND | 100 | 13 | ug/m3 |
| 1330-20-7 | 106.2 | Xylenes (total) | 30.4 | 24 | 3.0 | ppbv | | 132 | 100 | 13 | ug/m3 |

| CAS No. | Surrogate Recoveries | Run# 1 | Run# 2 | Limits |
|----------|----------------------|------------------|------------------|---------|
| 460-00-4 | 4-Bromofluorobenzene | 23% ^b | 22% ^b | 65-128% |

ND = Not detected MDL = Method Detection Limit
 RL = Reporting Limit
 E = Indicates value exceeds calibration range

J = Indicates an estimated value
 B = Indicates analyte found in associated method blank
 N = Indicates presumptive evidence of compound



Report of Analysis

Page 3 of 3

| | |
|--|--------------------------------|
| Client Sample ID: B8-SSV2/08312015 | Date Sampled: 08/31/15 |
| Lab Sample ID: JC2870-1A | Date Received: 09/02/15 |
| Matrix: AIR - Soil Vapor Comp. Summa ID: A1216,A681 | Percent Solids: n/a |
| Method: TO-15 | |
| Project: BMSMC, Building 5 Area, PR | |

VOA TO15 List

| CAS No. | MW | Compound | Result | RL | MDL | Units | Q | Result | RL | MDL | Units |
|---------|----|----------|--------|----|-----|-------|---|--------|----|-----|-------|
|---------|----|----------|--------|----|-----|-------|---|--------|----|-----|-------|

(a) Confirmation run.

(b) Outside control limits due to matrix interference.



ND = Not detected MDL = Method Detection Limit
RL = Reporting Limit
E = Indicates value exceeds calibration range

J = Indicates an estimated value
B = Indicates analyte found in associated method blank
N = Indicates presumptive evidence of a compound

Accutest LabLink@872636 09:45 14-Sep-2015

Report of Analysis

Page 1 of 3

| | | | |
|--------------------------|---------------------------------------|------------------------|----------|
| Client Sample ID: | B8-SSV1/08312015 | Date Sampled: | 08/31/15 |
| Lab Sample ID: | JC2870-2 | Date Received: | 09/02/15 |
| Matrix: | AIR - Soil Vapor Comp. Summa ID: A382 | Percent Solids: | n/a |
| Method: | TO-15 | | |
| Project: | BMSMC, Building 5 Area, PR | | |

| Run # | File ID | DF | Analyzed | By | Prep Date | Prep Batch | Analytical Batch |
|---------------------|-----------|----|----------|-----|-----------|------------|------------------|
| Run #1 | 3W49836.D | 1 | 09/03/15 | YMH | n/a | n/a | V3W1889 |
| Run #2 ^a | 3W49846.D | 1 | 09/03/15 | YMH | n/a | n/a | V3W1890 |


| | Initial Volume |
|--------|----------------|
| Run #1 | 100 ml |
| Run #2 | 100 ml |

VOA TO15 List

| CAS No. | MW | Compound | Result | RL | MDL | Units | Q | Result | RL | MDL | Units |
|------------|-------|----------------------------|--------|------|-------|-------|---|--------|-----|------|-------|
| 67-64-1 | 58.08 | Acetone | 106 | 0.80 | 0.13 | ppbv | | 252 | 1.9 | 0.31 | ug/m3 |
| 106-99-0 | 54.09 | 1,3-Butadiene | ND | 0.80 | 0.12 | ppbv | | ND | 1.8 | 0.27 | ug/m3 |
| 71-43-2 | 78.11 | Benzene | 1.5 | 0.80 | 0.12 | ppbv | | 4.8 | 2.6 | 0.38 | ug/m3 |
| 75-27-4 | 163.8 | Bromodichloromethane | ND | 0.80 | 0.13 | ppbv | | ND | 5.4 | 0.87 | ug/m3 |
| 75-25-2 | 252.8 | Bromoform | ND | 0.80 | 0.082 | ppbv | | ND | 8.3 | 0.85 | ug/m3 |
| 74-83-9 | 94.94 | Bromomethane | ND | 0.80 | 0.087 | ppbv | | ND | 3.1 | 0.34 | ug/m3 |
| 593-60-2 | 106.9 | Bromoethene | ND | 0.80 | 0.081 | ppbv | | ND | 3.5 | 0.35 | ug/m3 |
| 100-44-7 | 126 | Benzyl Chloride | ND | 0.80 | 0.10 | ppbv | | ND | 4.1 | 0.52 | ug/m3 |
| 75-15-0 | 76.14 | Carbon disulfide | ND | 0.80 | 0.11 | ppbv | | ND | 2.5 | 0.34 | ug/m3 |
| 108-90-7 | 112.6 | Chlorobenzene | ND | 0.80 | 0.13 | ppbv | | ND | 3.7 | 0.60 | ug/m3 |
| 75-00-3 | 64.52 | Chloroethane | ND | 0.80 | 0.086 | ppbv | | ND | 2.1 | 0.23 | ug/m3 |
| 67-66-3 | 119.4 | Chloroform | 0.69 | 0.80 | 0.12 | ppbv | J | 3.4 | 3.9 | 0.59 | ug/m3 |
| 74-87-3 | 50.49 | Chloromethane | 0.68 | 0.80 | 0.12 | ppbv | J | 1.4 | 1.7 | 0.25 | ug/m3 |
| 107-05-1 | 76.53 | 3-Chloropropene | ND | 0.80 | 0.11 | ppbv | | ND | 2.5 | 0.34 | ug/m3 |
| 95-49-8 | 126.6 | 2-Chlorotoluene | ND | 0.80 | 0.13 | ppbv | | ND | 4.1 | 0.67 | ug/m3 |
| 56-23-5 | 153.8 | Carbon tetrachloride | ND | 0.80 | 0.098 | ppbv | | ND | 5.0 | 0.62 | ug/m3 |
| 110-82-7 | 84.16 | Cyclohexane | 2.5 | 0.80 | 0.13 | ppbv | | 8.6 | 2.8 | 0.45 | ug/m3 |
| 75-34-3 | 98.96 | 1,1-Dichloroethane | ND | 0.80 | 0.12 | ppbv | | ND | 3.2 | 0.49 | ug/m3 |
| 75-35-4 | 96.94 | 1,1-Dichloroethylene | ND | 0.80 | 0.11 | ppbv | | ND | 3.2 | 0.44 | ug/m3 |
| 106-93-4 | 187.9 | 1,2-Dibromoethane | ND | 0.80 | 0.14 | ppbv | | ND | 6.1 | 1.1 | ug/m3 |
| 107-06-2 | 98.96 | 1,2-Dichloroethane | ND | 0.80 | 0.10 | ppbv | | ND | 3.2 | 0.40 | ug/m3 |
| 78-87-5 | 113 | 1,2-Dichloropropane | ND | 0.80 | 0.20 | ppbv | | ND | 3.7 | 0.92 | ug/m3 |
| 123-91-1 | 88.12 | 1,4-Dioxane | ND | 0.80 | 0.25 | ppbv | | ND | 2.9 | 0.90 | ug/m3 |
| 75-71-8 | 120.9 | Dichlorodifluoromethane | 0.48 | 0.80 | 0.15 | ppbv | J | 2.4 | 4.0 | 0.74 | ug/m3 |
| 124-48-1 | 208.3 | Dibromochloromethane | ND | 0.80 | 0.17 | ppbv | | ND | 6.8 | 1.4 | ug/m3 |
| 156-60-5 | 96.94 | trans-1,2-Dichloroethylene | ND | 0.80 | 0.081 | ppbv | | ND | 3.2 | 0.32 | ug/m3 |
| 156-59-2 | 96.94 | cis-1,2-Dichloroethylene | ND | 0.80 | 0.099 | ppbv | | ND | 3.2 | 0.39 | ug/m3 |
| 10061-01-5 | 111 | cis-1,3-Dichloropropene | ND | 0.80 | 0.14 | ppbv | | ND | 3.6 | 0.64 | ug/m3 |
| 541-73-1 | 147 | m-Dichlorobenzene | ND | 0.80 | 0.11 | ppbv | | ND | 4.8 | 0.66 | ug/m3 |
| 95-50-1 | 147 | o-Dichlorobenzene | ND | 0.80 | 0.12 | ppbv | | ND | 4.8 | 0.72 | ug/m3 |
| 106-46-7 | 147 | p-Dichlorobenzene | ND | 0.80 | 0.077 | ppbv | | ND | 4.8 | 0.46 | ug/m3 |
| 10061-02-6 | 111 | trans-1,3-Dichloropropene | ND | 0.80 | 0.082 | ppbv | | ND | 3.6 | 0.37 | ug/m3 |

ND = Not detected MDL = Method Detection Limit
 RL = Reporting Limit
 E = Indicates value exceeds calibration range

J = Indicates an estimated value
 B = Indicates analyte found in associated method blank
 N = Indicates presumptive evidence of a compound



Report of Analysis

Page 2 of 3

| | |
|--|--------------------------------|
| Client Sample ID: B8-SSV1/08312015 | Date Sampled: 08/31/15 |
| Lab Sample ID: JC2870-2 | Date Received: 09/02/15 |
| Matrix: AIR - Soil Vapor Comp. Summa ID: A382 | Percent Solids: n/a |
| Method: TO-15 | |
| Project: BMSMC, Building 5 Area, PR | |

4.3

4


VOA TO15 List

| CAS No. | MW | Compound | Result | RL | MDL | Units | Q | Result | RL | MDL | Units |
|-----------|--------|---------------------------|--------|------|-------|-------|---|--------|------|------|-------|
| 64-17-5 | 46.07 | Ethanol | 27.6 | 2.0 | 0.66 | ppbv | | 52.0 | 3.8 | 1.2 | ug/m3 |
| 100-41-4 | 106.2 | Ethylbenzene | 0.76 | 0.80 | 0.19 | ppbv | J | 3.3 | 3.5 | 0.83 | ug/m3 |
| 141-78-6 | 88 | Ethyl Acetate | 3.4 | 0.80 | 0.25 | ppbv | | 12 | 2.9 | 0.90 | ug/m3 |
| 622-96-8 | 120.2 | 4-Ethyltoluene | 0.36 | 0.80 | 0.088 | ppbv | J | 1.8 | 3.9 | 0.43 | ug/m3 |
| 76-13-1 | 187.4 | Freon 113 | ND | 0.80 | 0.11 | ppbv | | ND | 6.1 | 0.84 | ug/m3 |
| 76-14-2 | 170.9 | Freon 114 | ND | 0.80 | 0.10 | ppbv | | ND | 5.6 | 0.70 | ug/m3 |
| 142-82-5 | 100.2 | Heptane | ND | 0.80 | 0.12 | ppbv | | ND | 3.3 | 0.49 | ug/m3 |
| 87-68-3 | 260.8 | Hexachlorobutadiene | ND | 0.80 | 0.13 | ppbv | | ND | 8.5 | 1.4 | ug/m3 |
| 110-54-3 | 86.17 | Hexane | 0.86 | 0.80 | 0.11 | ppbv | | 3.0 | 2.8 | 0.39 | ug/m3 |
| 591-78-6 | 100 | 2-Hexanone | ND | 0.80 | 0.18 | ppbv | | ND | 3.3 | 0.74 | ug/m3 |
| 67-63-0 | 60.1 | Isopropyl Alcohol | 13.3 | 0.80 | 0.48 | ppbv | | 32.7 | 2.0 | 1.2 | ug/m3 |
| 75-09-2 | 84.94 | Methylene chloride | 0.80 | 0.80 | 0.54 | ppbv | | 2.8 | 2.8 | 1.9 | ug/m3 |
| 78-93-3 | 72.11 | Methyl ethyl ketone | 22.5 | 0.80 | 0.19 | ppbv | | 66.4 | 2.4 | 0.56 | ug/m3 |
| 108-10-1 | 100.2 | Methyl Isobutyl Ketone | ND | 0.80 | 0.11 | ppbv | | ND | 3.3 | 0.45 | ug/m3 |
| 1634-04-4 | 88.15 | Methyl Tert Butyl Ether | ND | 0.80 | 0.10 | ppbv | | ND | 2.9 | 0.36 | ug/m3 |
| 80-62-6 | 100.12 | Methylmethacrylate | ND | 0.80 | 0.12 | ppbv | | ND | 3.3 | 0.49 | ug/m3 |
| 115-07-1 | 42 | Propylene | 15.0 | 2.0 | 0.32 | ppbv | | 25.8 | 3.4 | 0.55 | ug/m3 |
| 100-42-5 | 104.1 | Styrene | ND | 0.80 | 0.10 | ppbv | | ND | 3.4 | 0.43 | ug/m3 |
| 71-55-6 | 133.4 | 1,1,1-Trichloroethane | ND | 0.80 | 0.13 | ppbv | | ND | 4.4 | 0.71 | ug/m3 |
| 79-34-5 | 167.9 | 1,1,2,2-Tetrachloroethane | ND | 0.80 | 0.12 | ppbv | | ND | 5.5 | 0.82 | ug/m3 |
| 79-00-5 | 133.4 | 1,1,2-Trichloroethane | ND | 0.80 | 0.14 | ppbv | | ND | 4.4 | 0.76 | ug/m3 |
| 120-82-1 | 181.5 | 1,2,4-Trichlorobenzene | ND | 0.80 | 0.18 | ppbv | | ND | 5.9 | 1.3 | ug/m3 |
| 95-63-6 | 120.2 | 1,2,4-Trimethylbenzene | 2.1 | 0.80 | 0.092 | ppbv | | 10 | 3.9 | 0.45 | ug/m3 |
| 108-67-8 | 120.2 | 1,3,5-Trimethylbenzene | 2.0 | 0.80 | 0.12 | ppbv | | 9.8 | 3.9 | 0.59 | ug/m3 |
| 540-84-1 | 114.2 | 2,2,4-Trimethylpentane | ND | 0.80 | 0.085 | ppbv | | ND | 3.7 | 0.40 | ug/m3 |
| 75-65-0 | 74.12 | Tertiary Butyl Alcohol | 1.1 | 0.80 | 0.20 | ppbv | | 3.3 | 2.4 | 0.61 | ug/m3 |
| 127-18-4 | 165.8 | Tetrachloroethylene | ND | 0.16 | 0.095 | ppbv | | ND | 1.1 | 0.64 | ug/m3 |
| 109-99-9 | 72.11 | Tetrahydrofuran | ND | 0.80 | 0.17 | ppbv | | ND | 2.4 | 0.50 | ug/m3 |
| 108-88-3 | 92.14 | Toluene | 1.4 | 0.80 | 0.081 | ppbv | | 5.3 | 3.0 | 0.31 | ug/m3 |
| 79-01-6 | 131.4 | Trichloroethylene | ND | 0.16 | 0.10 | ppbv | | ND | 0.86 | 0.54 | ug/m3 |
| 75-69-4 | 137.4 | Trichlorofluoromethane | ND | 0.80 | 0.081 | ppbv | | ND | 4.5 | 0.46 | ug/m3 |
| 75-01-4 | 62.5 | Vinyl chloride | ND | 0.80 | 0.13 | ppbv | | ND | 2.0 | 0.33 | ug/m3 |
| 108-05-4 | 86 | Vinyl Acetate | ND | 0.80 | 0.22 | ppbv | | ND | 2.8 | 0.77 | ug/m3 |
| | 106.2 | m,p-Xylene | 1.6 | 0.80 | 0.17 | ppbv | | 6.9 | 3.5 | 0.74 | ug/m3 |
| 95-47-6 | 106.2 | o-Xylene | 1.0 | 0.80 | 0.10 | ppbv | | 4.3 | 3.5 | 0.43 | ug/m3 |
| 1330-20-7 | 106.2 | Xylenes (total) | 2.6 | 0.80 | 0.10 | ppbv | | 11 | 3.5 | 0.43 | ug/m3 |

| CAS No. | Surrogate Recoveries | Run# 1 | Run# 2 | Limits |
|----------|----------------------|------------------|------------------|---------|
| 460-00-4 | 4-Bromofluorobenzene | 36% ^b | 39% ^b | 65-128% |

ND = Not detected MDL = Method Detection Limit
 RL = Reporting Limit
 E = Indicates value exceeds calibration range

J = Indicates an estimated value
 B = Indicates analyte found in associated method blank
 N = Indicates presumptive evidence of a compound



Report of Analysis

Page 3 of 3

| | |
|--|--------------------------------|
| Client Sample ID: B8-SSV1/08312015 | Date Sampled: 08/31/15 |
| Lab Sample ID: JC2870-2 | Date Received: 09/02/15 |
| Matrix: AIR - Soil Vapor Comp. Summa ID: A382 | Percent Solids: n/a |
| Method: TO-15 | |
| Project: BMSMC, Building 5 Area, PR | |

4.3

4

VOA TO15 List

| CAS No. | MW | Compound | Result | RL | MDL | Units | Q | Result | RL | MDL | Units |
|---------|----|----------|--------|----|-----|-------|---|--------|----|-----|-------|
|---------|----|----------|--------|----|-----|-------|---|--------|----|-----|-------|

(a) Confirmation run.

(b) Outside control limits due to matrix interference.



ND = Not detected MDL = Method Detection Limit
RL = Reporting Limit
E = Indicates value exceeds calibration range


J = Indicates an estimated value
B = Indicates analyte found in associated method blank
N = Indicates presumptive evidence of a compound

4

Page 1 of 1

| | |
|--------|-----------------------|
| | Initial Volume |
| Run #1 | 0.50 ml |
| Run #2 | |

| CAS No. | Surrogate Recoveries | Run# 1 | Run# 2 | Limits |
|----------|----------------------|--------|--------|---------|
| 460-00-4 | 4-Bromofluorobenzene | 95% | | 61-137% |
| 460-00-4 | 4-Bromofluorobenzene | 93% | | 61-137% |



J = Indicates an estimated value
B = Indicates analyte found in associated method blank
N = Indicates presumptive evidence of a compound

Accutest LabLink@872636 09:45 14-Sep-2015

Report of Analysis

Page 1 of 3

| | | | |
|-------------------|--|-----------------|----------|
| Client Sample ID: | B8-SSV1D/08312015 | Date Sampled: | 08/31/15 |
| Lab Sample ID: | JC2870-3 | Date Received: | 09/02/15 |
| Matrix: | AIR - Soil Vapor Comp. Summa ID: A1154 | Percent Solids: | n/a |
| Method: | TO-15 | | |
| Project: | BMSMC, Building 5 Area, PR | | |

| Run # | File ID | DF | Analyzed | By | Prep Date | Prep Batch | Analytical Batch |
|---------------------|-----------|----|----------|-----|-----------|------------|------------------|
| Run #1 | 3W49847.D | 1 | 09/03/15 | YMH | n/a | n/a | V3W1890 |
| Run #2 ^a | 3W49837.D | 1 | 09/03/15 | YMH | n/a | n/a | V3W1889 |


| | Initial Volume |
|--------|----------------|
| Run #1 | 50.0 ml |
| Run #2 | 20.0 ml |

VOA TO15 List

| CAS No. | MW | Compound | Result | RL | MDL | Units | Q | Result | RL | MDL | Units |
|------------|-------|----------------------------|--------|-----|------|-------|---|--------|------|-------|-------|
| 67-64-1 | 58.08 | Acetone | 88.3 | 1.6 | 0.25 | ppbv | | 210 | 3.8 | 0.59 | ug/m3 |
| 106-99-0 | 54.09 | 1,3-Butadiene | ND | 1.6 | 0.25 | ppbv | | ND | 3.5 | 0.55 | ug/m3 |
| 72-43-2 | 78.11 | Benzene | 1.4 | 1.6 | 0.24 | ppbv | J | 4.5 | 5.1 | 0.77 | ug/m3 |
| 75-27-4 | 163.8 | Bromodichloromethane | ND | 1.6 | 0.25 | ppbv | | ND | 11 | 1.7 | ug/m3 |
| 75-25-2 | 252.8 | Bromoform | ND | 1.6 | 0.16 | ppbv | | ND | 17 | 1.7 | ug/m3 |
| 74-83-9 | 94.94 | Bromomethane | ND | 1.6 | 0.17 | ppbv | | ND | 6.2 | 0.66 | ug/m3 |
| 593-60-2 | 106.9 | Bromoethene | ND | 1.6 | 0.16 | ppbv | | ND | 7.0 | 0.70 | ug/m3 |
| 100-44-7 | 126 | Benzyl Chloride | ND | 1.6 | 0.21 | ppbv | | ND | 8.2 | 1.1 | ug/m3 |
| 75-15-0 | 76.14 | Carbon disulfide | ND | 1.6 | 0.23 | ppbv | | ND | 5.0 | 0.72 | ug/m3 |
| 108-90-7 | 112.6 | Chlorobenzene | ND | 1.6 | 0.26 | ppbv | | ND | 7.4 | 1.2 | ug/m3 |
| 75-00-3 | 64.52 | Chloroethane | ND | 1.6 | 0.17 | ppbv | | ND | 4.2 | 0.45 | ug/m3 |
| 67-66-3 | 119.4 | Chloroform | ND | 1.6 | 0.25 | ppbv | | ND | 7.8 | 1.2 | ug/m3 |
| 74-87-3 | 50.49 | Chloromethane | ND | 1.6 | 0.23 | ppbv | | ND | 3.3 | 0.47 | ug/m3 |
| 107-05-1 | 76.53 | 3-Chloropropene | ND | 1.6 | 0.22 | ppbv | | ND | 5.0 | 0.69 | ug/m3 |
| 95-49-8 | 126.6 | 2-Chlorotoluene | ND | 1.6 | 0.26 | ppbv | | ND | 8.3 | 1.3 | ug/m3 |
| 56-23-5 | 153.8 | Carbon tetrachloride | ND | 1.6 | 0.20 | ppbv | | ND | 10 | 1.3 | ug/m3 |
| 110-82-7 | 84.16 | Cyclohexane | 2.2 | 1.6 | 0.26 | ppbv | | 7.6 | 5.5 | 0.89 | ug/m3 |
| 75-34-3 | 98.96 | 1,1-Dichloroethane | ND | 1.6 | 0.25 | ppbv | | ND | 6.5 | 1.0 | ug/m3 |
| 75-35-4 | 96.94 | 1,1-Dichloroethylene | ND | 1.6 | 0.23 | ppbv | | ND | 6.3 | 0.91 | ug/m3 |
| 106-93-4 | 187.9 | 1,2-Dibromoethane | ND | 1.6 | 0.28 | ppbv | | ND | 12 | 2.2 | ug/m3 |
| 107-06-2 | 98.96 | 1,2-Dichloroethane | ND | 1.6 | 0.21 | ppbv | | ND | 6.5 | 0.85 | ug/m3 |
| 78-87-5 | 113 | 1,2-Dichloropropane | ND | 1.6 | 0.40 | ppbv | | ND | 7.4 | 1.8 | ug/m3 |
| 123-91-1 | 88.12 | 1,4-Dioxane | ND | 1.6 | 0.50 | ppbv | | ND | 5.8 | 1.8 | ug/m3 |
| 75-71-8 | 120.9 | Dichlorodifluoromethane | ND | 1.6 | 0.29 | ppbv | | ND | 7.9 | 1.4 | ug/m3 |
| 124-48-1 | 208.3 | Dibromochloromethane | ND | 1.6 | 0.33 | ppbv | | ND | 14 | 2.8 | ug/m3 |
| 156-60-5 | 96.94 | trans-1,2-Dichloroethylene | ND | 1.6 | 0.16 | ppbv | | ND | 6.3 | 0.63 | ug/m3 |
| 156-59-2 | 96.94 | cis-1,2-Dichloroethylene | ND | 1.6 | 0.20 | ppbv | | ND | 6.3 | 0.79 | ug/m3 |
| 10061-01-5 | 111 | cis-1,3-Dichloropropene | ND | 1.6 | 0.28 | ppbv | | ND | 7.3 | 1.3 | ug/m3 |
| 541-73-1 | 147 | m-Dichlorobenzene | ND | 1.6 | 0.22 | ppbv | | ND | 9.6 | 1.3 | ug/m3 |
| 95-50-1 | 147 | o-Dichlorobenzene | ND | 1.6 | 0.24 | ppbv | | ND | 1.4 | ug/m3 | |
| 106-46-7 | 147 | p-Dichlorobenzene | ND | 1.6 | 0.15 | ppbv | | ND | 0.90 | ug/m3 | |
| 10061-02-6 | 111 | trans-1,3-Dichloropropene | ND | 1.6 | 0.16 | ppbv | | ND | 7.3 | 0.73 | ug/m3 |

ND = Not detected MDL = Method Detection Limit
 RL = Reporting Limit
 E = Indicates value exceeds calibration range

J = Indicates an estimated value
 B = Indicates analyte found in associated method blank
 N = Indicates presumptive evidence of compound



Report of Analysis

Page 2 of 3

| | | | |
|--------------------------|--|------------------------|----------|
| Client Sample ID: | B8-SSV1D/08312015 | Date Sampled: | 08/31/15 |
| Lab Sample ID: | JC2870-3 | Date Received: | 09/02/15 |
| Matrix: | AIR - Soil Vapor Comp. Summa ID: A1154 | Percent Solids: | n/a |
| Method: | TO-15 | | |
| Project: | BMSMC, Building 5 Area, PR | | |


VOA TO15 List

| CAS No. | MW | Compound | Result | RL | MDL | Units | Q | Result | RL | MDL | Units |
|-----------|--------|---------------------------|--------|------|------|-------|---|--------|-----|------|-------|
| 64-17-5 | 46.07 | Ethanol | 24.1 | 4.0 | 1.3 | ppbv | | 45.4 | 7.5 | 2.4 | ug/m3 |
| 100-41-4 | 106.2 | Ethylbenzene | 1.7 | 1.6 | 0.38 | ppbv | | 7.4 | 6.9 | 1.7 | ug/m3 |
| 141-78-6 | 88 | Ethyl Acetate | 3.0 | 1.6 | 0.51 | ppbv | | 11 | 5.8 | 1.8 | ug/m3 |
| 622-96-8 | 120.2 | 4-Ethyltoluene | 0.87 | 1.6 | 0.18 | ppbv | J | 4.3 | 7.9 | 0.88 | ug/m3 |
| 76-13-1 | 187.4 | Freon 113 | ND | 1.6 | 0.22 | ppbv | | ND | 12 | 1.7 | ug/m3 |
| 76-14-2 | 170.9 | Freon 114 | ND | 1.6 | 0.20 | ppbv | | ND | 11 | 1.4 | ug/m3 |
| 142-82-5 | 100.2 | Heptane | ND | 1.6 | 0.24 | ppbv | | ND | 6.6 | 0.98 | ug/m3 |
| 87-68-3 | 260.8 | Hexachlorobutadiene | ND | 1.6 | 0.26 | ppbv | | ND | 17 | 2.8 | ug/m3 |
| 110-54-3 | 86.17 | Hexane | ND | 1.6 | 0.22 | ppbv | | ND | 5.6 | 0.78 | ug/m3 |
| 591-78-6 | 100 | 2-Hexanone | ND | 1.6 | 0.35 | ppbv | | ND | 6.5 | 1.4 | ug/m3 |
| 67-63-0 | 60.1 | Isopropyl Alcohol | 12.7 | 1.6 | 0.96 | ppbv | | 31.2 | 3.9 | 2.4 | ug/m3 |
| 75-09-2 | 84.94 | Methylene chloride | 1.6 | 1.6 | 1.1 | ppbv | | 5.6 | 5.6 | 3.8 | ug/m3 |
| 78-93-3 | 72.11 | Methyl ethyl ketone | 18.0 | 1.6 | 0.39 | ppbv | | 53.1 | 4.7 | 1.2 | ug/m3 |
| 108-10-1 | 100.2 | Methyl Isobutyl Ketone | ND | 1.6 | 0.22 | ppbv | | ND | 6.6 | 0.90 | ug/m3 |
| 1634-04-4 | 88.15 | Methyl Tert Butyl Ether | ND | 1.6 | 0.21 | ppbv | | ND | 5.8 | 0.76 | ug/m3 |
| 80-62-6 | 100.12 | Methylmethacrylate | ND | 1.6 | 0.24 | ppbv | | ND | 6.6 | 0.98 | ug/m3 |
| 115-07-1 | 42 | Propylene | 20.9 | 4.0 | 0.65 | ppbv | | 35.9 | 6.9 | 1.1 | ug/m3 |
| 100-42-5 | 104.1 | Styrene | ND | 1.6 | 0.21 | ppbv | | ND | 6.8 | 0.89 | ug/m3 |
| 71-55-6 | 133.4 | 1,1,1-Trichloroethane | ND | 1.6 | 0.25 | ppbv | | ND | 8.7 | 1.4 | ug/m3 |
| 79-34-5 | 167.9 | 1,1,2,2-Tetrachloroethane | ND | 1.6 | 0.24 | ppbv | | ND | 11 | 1.6 | ug/m3 |
| 79-00-5 | 133.4 | 1,1,2-Trichloroethane | ND | 1.6 | 0.28 | ppbv | | ND | 8.7 | 1.5 | ug/m3 |
| 120-82-1 | 181.5 | 1,2,4-Trichlorobenzene | ND | 1.6 | 0.35 | ppbv | | ND | 12 | 2.6 | ug/m3 |
| 95-63-6 | 120.2 | 1,2,4-Trimethylbenzene | 3.7 | 1.6 | 0.18 | ppbv | | 18 | 7.9 | 0.88 | ug/m3 |
| 108-67-8 | 120.2 | 1,3,5-Trimethylbenzene | 3.9 | 1.6 | 0.24 | ppbv | | 19 | 7.9 | 1.2 | ug/m3 |
| 540-84-1 | 114.2 | 2,2,4-Trimethylpentane | ND | 1.6 | 0.17 | ppbv | | ND | 7.5 | 0.79 | ug/m3 |
| 75-65-0 | 74.12 | Tertiary Butyl Alcohol | 1.4 | 1.6 | 0.40 | ppbv | J | 4.2 | 4.9 | 1.2 | ug/m3 |
| 127-18-4 | 165.8 | Tetrachloroethylene | ND | 0.32 | 0.19 | ppbv | | ND | 2.2 | 1.3 | ug/m3 |
| 109-99-9 | 72.11 | Tetrahydrofuran | ND | 1.6 | 0.35 | ppbv | | ND | 4.7 | 1.0 | ug/m3 |
| 108-88-3 | 92.14 | Toluene | 2.2 | 1.6 | 0.16 | ppbv | | 8.3 | 6.0 | 0.60 | ug/m3 |
| 79-01-6 | 131.4 | Trichloroethylene | ND | 0.32 | 0.20 | ppbv | | ND | 1.7 | 1.1 | ug/m3 |
| 75-69-4 | 137.4 | Trichlorofluoromethane | ND | 1.6 | 0.16 | ppbv | | ND | 9.0 | 0.90 | ug/m3 |
| 75-01-4 | 62.5 | Vinyl chloride | ND | 1.6 | 0.26 | ppbv | | ND | 4.1 | 0.66 | ug/m3 |
| 108-05-4 | 86 | Vinyl Acetate | ND | 1.6 | 0.44 | ppbv | | ND | 5.6 | 1.5 | ug/m3 |
| | 106.2 | m,p-Xylene | 3.8 | 1.6 | 0.35 | ppbv | | 17 | 6.9 | 1.5 | ug/m3 |
| 95-47-6 | 106.2 | o-Xylene | 1.7 | 1.6 | 0.21 | ppbv | | 7.4 | 6.9 | 0.91 | ug/m3 |
| 1330-20-7 | 106.2 | Xylenes (total) | 5.5 | 1.6 | 0.21 | ppbv | | 24 | 6.9 | 0.91 | ug/m3 |

| CAS No. | Surrogate Recoveries | Run# 1 | Run# 2 | Limits |
|----------|----------------------|------------------|--------|---------|
| 460-00-4 | 4-Bromofluorobenzene | 55% ^b | 82% | 65-128% |

ND = Not detected MDL = Method Detection Limit
 RL = Reporting Limit
 E = Indicates value exceeds calibration range

J = Indicates an estimated value
 B = Indicates analyte found in associated method blank
 N = Indicates presumptive evidence of compound



Report of Analysis

Page 3 of 3

| | |
|---|--------------------------------|
| Client Sample ID: B8-SSV1D/08312015 | Date Sampled: 08/31/15 |
| Lab Sample ID: JC2870-3 | Date Received: 09/02/15 |
| Matrix: AIR - Soil Vapor Comp. Summa ID: A1154 | Percent Solids: n/a |
| Method: TO-15 | |
| Project: BMSMC, Building 5 Area, PR | |

4.4

4

VOA TO15 List

| CAS No. | MW | Compound | Result | RL | MDL | Units | Q | Result | RL | MDL | Units |
|---------|----|----------|--------|----|-----|-------|---|--------|----|-----|-------|
|---------|----|----------|--------|----|-----|-------|---|--------|----|-----|-------|

(a) Confirmation run.

(b) Outside control limits due to matrix interference.



ND = Not detected MDL = Method Detection Limit
RL = Reporting Limit
E = Indicates value exceeds calibration range

J = Indicates an estimated value
B = Indicates analyte found in associated method blank
N = Indicates presumptive evidence of a compound

Accutest LabLink@872636 09:45 14-Sep-2015

Report of Analysis

Page 1 of 1

| | | | |
|--------------------------|--|------------------------|----------|
| Client Sample ID: | B8-SSV1D/08312015 | Date Sampled: | 08/31/15 |
| Lab Sample ID: | JC2870-3 | Date Received: | 09/02/15 |
| Matrix: | AIR - Soil Vapor Comp. Summa ID: A1154 | Percent Solids: | n/a |
| Method: | EPA TO-3 | | |
| Project: | BMSMC, Building 5 Area, PR | | |

| | File ID | DF | Analyzed | By | Prep Date | Prep Batch | Analytical Batch |
|--------|-----------|------|----------|-----|-----------|------------|------------------|
| Run #1 | QT19209.D | 1.52 | 09/04/15 | TCH | n/a | n/a | GQT1157 |
| Run #2 | | | | | | | |

| | Initial Volume |
|--------|----------------|
| Run #1 | 0.50 ml |
| Run #2 | |

| CAS No. | MW | Compound | Result | RL | MDL | Units | Q | Result | RL | MDL | Units |
|----------|----|----------|--------|------|-------|-------|---|--------|-----|------|-------|
| 74-82-8 | 16 | Methane | 197 | 7.6 | 0.51 | ppmv | | 129 | 5.0 | 0.33 | mg/m3 |
| 74-98-6 | 44 | Propane | ND | 0.76 | 0.063 | ppmv | | ND | 1.4 | 0.11 | mg/m3 |
| 106-97-8 | 58 | n-Butane | ND | 0.76 | 0.091 | ppmv | | ND | 1.8 | 0.22 | mg/m3 |

| CAS No. | Surrogate Recoveries | Run# 1 | Run# 2 | Limits |
|----------|----------------------|--------|--------|---------|
| 460-00-4 | 4-Bromofluorobenzene | 102% | | 61-137% |
| 460-00-4 | 4-Bromofluorobenzene | 103% | | 61-137% |



ND = Not detected MDL = Method Detection Limit
 RL = Reporting Limit
 E = Indicates value exceeds calibration range

J = Indicates an estimated value
 B = Indicates analyte found in associated method blank
 N = Indicates presumptive evidence of a compound

CHAIN OF CUSTODY

Air Sampling Field Data Sheet

801279535066

MS-8/18/2015-7

PAGE 1 OF 1

Lab Job # 7E2870

[illegible]

5

JC2870: Chain of Custody

Page 1 of 3

DATA REVIEW WORKSHEETS

Project Number: jc2870-1

Date: 08/31/2015

REVIEW OF VOLATILE ORGANIC PACKAGE

The following guidelines for evaluating volatile organics were created to delineate required validation actions. This document will assist the reviewer in using professional judgment to make more informed decision and in better serving the needs of the data users. The sample results were assessed according to USEPA data validation guidance documents in the following order of precedence: QC criteria from "Compendium Method TO-15. Determination of Volatile Organic Compounds (VOCs) In Air Collected In Specially-Prepared Canisters and Analyzed By Gas Chromatography/Mass Spectrometry (GC/MS), January, 1999"; USEPA Hazardous Waste Support Branch. Validating Air Samples. Volatile Organic Analysis of Ambient Air in Canisters by Method TO-15, (SOP # HW-31. Revision #4. October, 2006). The QC criteria and data validation actions listed on the data review worksheets are from the primary guidance document, unless otherwise noted.

The hardcopied (laboratory name) Eurofins - Air Toxics data package received has been reviewed and the quality control and performance data summarized. The data review for VOCs included:

Lab. Project/SDG No.: jc2870-1

Sample matrix: Air

No. of Samples: 3

Trip blank No.: -

Field blank No.: -

Equipment blank No.: -

Field duplicate No.: B8-SSV1/B8-SSV2

☒ Data Completeness

☒ Laboratory Control Spikes

☒ Holding Times

☒ Field Duplicates

☒ GC/MS Tuning

☒ Calibrations

☒ Internal Standard Performance

☒ Compound Identifications

☒ Blanks

☒ Compound Quantitation

☒ Surrogate Recoveries

☒ Quantitation Limits

☐ N/A Matrix Spike/Matrix Spike Duplicate

Overall Comments: VOCs by method TO-15

Definition of Qualifiers:

J- Estimated results

U- Compound not detected

R- Rejected data

UJ- Estimated nondetect

Reviewer: Rafael Infante

Date: 09/15/2015

DATA COMPLETENESS

DATE RECEIVED

A thick, solid gray diagonal line runs from the top-left corner of the page towards the bottom-right corner. It is positioned such that it divides the page into two roughly equal triangular sections. The line is uniform in thickness and has sharp, clean edges.

DATA REVIEW WORKSHEETS

All criteria were met X
Criteria were not met
and/or see below

HOLDING TIMES

The objective of this parameter is to ascertain the validity of the results based on the holding time of the sample from time of collection to the time of analysis.

Complete table for all samples and note the analysis and/or preservation not within criteria

| SAMPLE ID | DATE SAMPLED | DATE ANALYZED | pH | ACTION |
|---|--------------|---------------|----|--------|
| | | | | |
| | | | | |
| All samples analyzed within the recommended method holding time | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |

Criteria

Aqueous samples – 14 days from sample collection for preserved samples ($\text{pH} \leq 2$, 4°C), no air bubbles.

Aqueous samples – 7 days from sample collection for unpreserved samples, 4°C , no air bubbles.

Soil samples- 7 days from sample collection.

Cooler temperature (Criteria: $4 \pm 2^{\circ}\text{C}$): N/A – summa canisters

Actions

If the VOCs vial(s) have air bubbles, estimate positive results (J) and reject nondetects (R).

If the % solids of soil samples is 10-50%, estimate positive results (J) and nondetects (UJ).

If the % solid of soil samples is $< 10\%$, estimate positive results (J) and reject nondetects (R).

If holding times are exceeded but < 14 days beyond criteria, estimate positive results (J) and nondetects (UJ).

If holding times are exceeded but < 28 days beyond criteria, estimate positive results (J) and reject nondetects (R).

If holding times are grossly exceeded (> 28 days beyond criteria), reject all results (R).

If samples were not iced or if the ice were melted ($> 10^{\circ}\text{C}$), estimate positive results (J) and nondetects (UJ).

GC/MS TUNING

The assessment of the tuning results is to determine if the sample instrumentation is within the standard tuning QC limits

 X BFB tuning was performed for every 24 hours of sample analysis.

If no, use professional judgment to determine whether the associated data should be accepted, qualified or rejected.

List the samples affected:

If mass calibration is in error, all associated data are rejected.

DATA REVIEW WORKSHEETS

All criteria were met X
 Criteria were not met
 and/or see below

CALIBRATION VERIFICATION

Compliance requirements for satisfactory instrument calibration are established to ensure that the instrument is capable of producing and maintaining acceptable quantitative data.

Date of initial calibration: 07/16/15; 08/14/15; 08/28/15
 Dates of continuing calibration: 07/29/15; 08/29/15; 09/02/15 09/03/15
 Instrument ID numbers: GCMS3W
 Matrix/Level: Air/low

| DATE | LAB ID# | FILE | CRITERIA OUT RFs, %RSD, %D, r | COMPOUND | SAMPLES AFFECTED |
|---|------------|------|----------------------------------|----------|---------------------|
| | | | | | |
| | | | | | |
| | | | | | |
| Initial and continuing calibration met the method performance criteria. | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |

Criteria

All RFs must be > 0.05 regardless of method requirements for SPCC.

All %RSD must be $\leq 15\%$ regardless of method requirements for CCC.

All %Ds must be $\leq 30\%$ regardless of method requirements for CCC.

Method TO-15 does not specify criterion for the curve correlation coefficient (r). A limit for r of ≥ 0.995 has therefore been utilized as professional judgment.

Actions

If any compound has an initial RF or a continuing RF of < 0.05 , estimate positive results (J) and reject nondetects (R), regardless of method requirements.

If any compound has a %RSD $> 15\%$, estimate positive results (J) and use professional judgment to qualify nondetects.

If any compound has a %RSD $> 90\%$, estimate positive results (J) and reject nondetects (R).

If any compound has a % D $> 30\%$, estimate positive results (J) and reject nondetects (R).

If any compound has a % D $> 30\%$, estimate positive results (J) and nondetects (UJ).

If any compound has a % D $> 90\%$, estimate positive results (J) and reject nondetects (R).

If any compound has $r < 0.995$, estimate positive results and nondetects.

A separate worksheet should be filled for each initial curve

DATA REVIEW WORKSHEETS

All criteria were met X
 Criteria were not met
 and/or see below _____

CALIBRATION VERIFICATION

Compliance requirements for satisfactory instrument calibration are established to ensure that the instrument is capable of producing and maintaining acceptable quantitative data.

Date of initial calibration: 08/06/15
 Dates of continuing calibration: 08/07/15
 Instrument ID numbers: GCMS5W
 Matrix/Level: Air/low

| DATE | LAB ID# | FILE | CRITERIA OUT RFs, %RSD, %D, r | COMPOUND | SAMPLES AFFECTED |
|---|------------|------|----------------------------------|----------|---------------------|
| | | | | | |
| | | | | | |
| | | | | | |
| Initial and continuing calibration met the method performance criteria. | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |

Criteria

All RFs must be > 0.05 regardless of method requirements for SPCC.
 All %RSD must be $\leq 15\%$ regardless of method requirements for CCC.
 All %Ds must be $\leq 30\%$ regardless of method requirements for CCC.
 Method TO-15 does not specify criterion for the curve correlation coefficient (r). A limit for r of ≥ 0.995 has therefore been utilized as professional judgment.

Actions

If any compound has an initial RF or a continuing RF of < 0.05 , estimate positive results (J) and reject nondetects (R), regardless of method requirements.
 If any compound has a %RSD $> 15\%$, estimate positive results (J) and use professional judgment to qualify nondetects.
 If any compound has a %RSD $> 90\%$, estimate positive results (J) and reject nondetects (R).
 If any compound has a % D $> 30\%$, estimate positive results (J) and reject nondetects (R).
 If any compound has a % D $> 30\%$, estimate positive results (J) and nondetects (UJ).
 If any compound has a % D $> 90\%$, estimate positive results (J) and reject nondetects (R).
 If any compound has $r < 0.995$, estimate positive results and nondetects.

A separate worksheet should be filled for each initial curve

DATA REVIEW WORKSHEETS

All criteria were met X
Criteria were not met
and/or see below _____

V A. BLANK ANALYSIS RESULTS (Sections 1 & 2)

The assessment of the blank analysis results is to determine the existence and magnitude of contamination problems. The criteria for evaluation of blanks apply only to blanks associated with the samples, including trip, equipment, and laboratory blanks. If problems with any blanks exist, all data associated with the case must be carefully evaluated to determine whether or not there is an inherent variability in the data for the case, or if the problem is an isolated occurrence not affecting other data.

List the contamination in the blanks below. High and low levels blanks must be treated separately.

Laboratory blanks

| DATE ANALYZED | LAB ID | LEVEL/ MATRIX | COMPOUND | CONCENTRATION/ UNITS |
|---------------|--------|------------------|----------|-------------------------|
|---------------|--------|------------------|----------|-------------------------|

All_method_blank_meeth_method_specific_criteria

Summa canisters met cleaning certification criteria

Field/Equipment/Trip blank

| DATE ANALYZED | LAB ID | LEVEL/ MATRIX | COMPOUND | CONCENTRATION UNITS |
|---------------|--------|------------------|----------|------------------------|
|---------------|--------|------------------|----------|------------------------|

No field/trip/equipment blanks analyzed with this data package.

DATA REVIEW WORKSHEETS

All criteria were met X
 Criteria were not met
 and/or see below

V B. BLANK ANALYSIS RESULTS (Section 3)

Blank Actions

Action Levels (ALs) should be based upon the highest concentration of contaminant determined in any blank. Do not qualify any blank with another blank. The ALs for samples which have been diluted should be corrected for the sample dilution factor and/or % moisture, where applicable. No positive sample results should be reported unless the concentration of the compound in the samples exceeds the ALs:

ALs = 10x the amount of common contaminants (methylene chloride, acetone, 2-butanone, and toluene)

ALs = 5x for any other compounds

Specific actions are as follows:

If the concentration is < sample quantitation limit (SQL) and \leq AL, report the compound as not detected (U) at the SQL.

If the concentration is \geq SQL but \leq AL, report the compound as not detected (U) at the reported concentration.

If the concentration is \geq SQL and $>$ AL, report the concentration unqualified.

Notes:

High and low level blanks must be treated separately

Compounds qualified "U" for blank contamination are still considered "hits" when qualifying for calibration criteria.

| CONTAMINATION SOURCE/LEVEL | COMPOUND | CONC/UNITS | AL/UNITS | SQL | AFFECTED SAMPLES |
|-------------------------------|----------|------------|----------|-----|---------------------|
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |

DATA REVIEW WORKSHEETS

All criteria were met _____
 Criteria were not met _____
 and/or see below _____X_____

SURROGATE SPIKE RECOVERIES

Laboratory performance of individual samples is established by evaluation of surrogate spike recoveries. All samples are spiked with surrogate compounds prior to sample analysis. The accuracy of the analysis is measured by the surrogate percent recovery. Since the effects of the sample matrix are frequently outside the control of the laboratory and may present relatively unique problems, the validation of data is frequently subjective and demands analytical experience and professional judgment.

List the percent recoveries (%Rs) which do not meet the criteria for surrogate recovery.

Matrix: solid/aqueous

SAMPLE ID SURROGATE COMPOUND ACTION

1,2-DICHLOROETHANE- Toluene- 4-BFB
 d4 d8

4-BFB used as surrogate; recoveries within laboratory control limits except for the followings: _____

| | | |
|-----------|------|-----------|
| JC2870-2 | 39 % | No action |
| JC2870-2 | 36 % | No action |
| JC2870-3 | 55 % | No action |
| JC2870-1A | 22 % | No action |
| JC2870-1A | 23 % | No action |

QC Limits* (Air)

LL to UL to to 65 to 128

- * QC limits are laboratory in-house performance criteria, LL = lower limit, UL = upper limit.
- * If QC limits are not available, use limits of 80 – 120 % for aqueous and 70 – 130 % for solid samples.

Actions:

| QUALITY | %R < 10% | %R = 10% - LL | %R > UL |
|--------------------|----------|---------------|---------|
| Positive results | J | J | J |
| Nondetects results | R | UJ | Accept |

Surrogate action should be applied:

If one or more surrogate in the VOC fraction is out of specification, but has a recovery of > 10%.

If any one surrogate in a fraction shows < 10 % recovery.

DATA REVIEW WORKSHEETS

All criteria were met _____
 Criteria were not met _____
 and/or see below __N/A__

VII. A MATRIX SPIKE/MATRIX SPIKE DUPLICATE (MS/MSD)

This data is generated to determine long term precision and accuracy in the analytical method for various matrices. This data alone cannot be used to evaluate the precision and accuracy of individual samples. If any % R in the MS or MSD falls outside the designated range, the reviewer should determine if there are matrix effects, i.e. LCS data are within the QC limits but MS/MSD data are outside QC limit.

1. MS/MSD Recoveries and Precision Criteria

The laboratory should use one MS and a duplicate analysis of an unspiked field sample if target analytes are expected in the sample. If target analytes are not expected, MS/MSD should be analyzed.

List the %Rs, RPD of the compounds which do not meet the criteria.

Sample ID: _____ Matrix/Level: _____

| MS OR MSD | COMPOUND | % R | RPD | QC LIMITS | ACTION |
|--|----------|-----|-----|-----------|--------|
| _____ MS/MSD are not required as part of Method TO-15; blank spike used to assess _____ accuracy _____ | | | | | |

- * QC limits are laboratory in-house performance criteria, LL = lower limit, UL = upper limit.
- * If QC limits are not available, use limits of 70 – 130 %.

Actions:

| QUALITY | %R < LL | %R > UL |
|--------------------|---------|---------|
| Positive results | J | J |
| Nondetects results | R | Accept |

MS/MSD criteria apply only to the unspiked sample, its dilutions, and the associated MS/MSD samples:

If the % R for the affected compounds were < LL (or 70 %), qualify positive results (J) and nondetects (UJ).

If the % R for the affected compounds were > UL (or 130 %), only qualify positive results (J).

If 25 % or more of all MS/MSD %R were < LL (or 70 %) or if two or more MS/MSD %Rs were < 10%, qualify all positive results (J) and reject nondetects (R).

A separate worksheet should be used for each MS/MSD pair.

DATA REVIEW WORKSHEETS

All criteria were met _____
Criteria were not met _____
and/or see below _____ N/A _____

VII. B MATRIX SPIKE/MATRIX SPIKE DUPLICATE

MS/MSD – Unspiked Compounds

It should be noted that Method TO-15 does not specify a MS/MSD criteria for the unspiked compounds in the sample. A %RSD of < 50% has therefore been utilized as professional judgment.

If all target analytes were spiked in the MS/MSD, this review element is not applicable.

List the %RSD of the compounds which do not meet the criteria.

Sample ID: _____ Matrix/Level/Unit: _____

[illegible]

Actions:

- * If the % RSD > 50, qualify the positive result in the unspiked samples as estimated (J).
* If the % RSD is not calculated (NC) due to nondetected value, use professional judgment to qualify the data.

DATA REVIEW WORKSHEETS

All criteria were met _____
 Criteria were not met _____
 and/or see below ____X____

VIII. LABORATORY CONTROL SAMPLE (LCS) ANALYSIS

This data is generated to determine accuracy of the analytical method for various matrices.

1. LCS Recoveries Criteria

Where LCS spiked with the same analyte at the same concentrations as the MS/MSD?
 Yes or No. If no make note in data review memo.

List the %R of compounds which do not meet the criteria

| LCS ID | COMPOUND | % R | QC LIMIT |
|--|----------|-----|----------|
| Blank_spike/spike_duplicate_%_recoveries_and_RPD_within_laboratory_control_limits_____ | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |

- * QC limits are laboratory in-house performance criteria, LL = lower limit, UL = upper limit.
- * If QC limits are not available, use limits of 70 – 130 %.

Actions:

| QUALITY | %R < LL | %R > UL |
|--------------------|---------|---------|
| Positive results | J | J |
| Nondetects results | R | Accept |

All analytes in the associated sample results are qualified for the following criteria.

If 25 % of the LCS recoveries were < LL (or 70 %), qualify all positive results (j) and reject nondetects (R).
 If two or more LCS were below 10 %, qualify all positive results as (J) and reject nondetects (R).

2. Frequency Criteria:

Where LCS analyzed at the required frequency and for each matrix? Yes or No.

If no, the data may be affected. Use professional judgment to determine the severity of the effect and qualify data accordingly. Discuss any actions below and list the samples affected.

DATA REVIEW WORKSHEETS

All criteria were met X
 Criteria were not met
 and/or see below

IX. LABORATORY DUPLICATE PRECISION

Sample IDs: JC2480-6/-6_Duplicate
 Sample IDs: JC1343-1/-1_Duplicate

Matrix: Air
 Matrix: Air

Laboratory duplicates samples may be taken and analyzed as an indication of overall precision. These analyses measure both field and lab precision; therefore, the results may have more variability than laboratory duplicates which only laboratory performance. It is also expected that soil duplicate results will have a greater variance than water matrices due to difficulties associated with collecting identical field duplicate samples.

The project QAPP should be reviewed for project-specific information.

Suggested criteria: RPD \pm 25% for air samples. If both samples and duplicate are <5 SQL, the RPD criteria is doubled.

| COMPOUND | SQL | SAMPLE CONC. | DUPLICATE CONC. | RPD | ACTION |
|--|-----|--------------|-----------------|-----|--------|
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| RPD within method performance criteria for all analytes. | | | | | |
| | | | | | |

Actions:

Qualify as estimated positive results (J) and nondetects (UJ) for the compound that exceeded the above criteria. For organics, only the sample and duplicate will be qualified.

If an RPD cannot be calculated because one or both of the sample results is not detected, the following actions apply:

If one sample result is not detected and the other is greater than 5x the SQL qualify (J/UJ).

If one sample value is not detected and the other is greater than 5x the SQL and the SQLs for the sample and duplicate are significantly different, use professional judgment to determine if qualification is appropriate.

If one sample value is not detected and the other is less than 5x, use professional judgment to determine if qualification is appropriate.

If both sample and duplicate results are not detected, no action is needed.

DATA REVIEW WORKSHEETS

All criteria were met _____
 Criteria were not met _____
 and/or see below X

IX. FIELD DUPLICATE PRECISION

Sample IDs: B8SSV1/B8SSV1D

Matrix: Air

Field duplicate samples may be taken and analyzed as an indication of overall precision. These analyses measure both field and lab precision; therefore, the results may have more variability than laboratory duplicates which only laboratory performance. It is also expected that soil duplicate results will have a greater variance than water matrices due to difficulties associated with collecting identical field duplicate samples.

The project QAPP should be reviewed for project-specific information.

Suggested criteria: RPD \pm 25% for air samples. If both samples and duplicate are <5 SQL, the RPD criteria is doubled.

| COMPOUND | RL | SAMPLE CONC. | DUPLICATE CONC. | RPD | ACTION |
|---|----|--------------|-----------------|-----|--------|
| | | | | | |
| | | | | | |
| RPD WITHIN GENERALLY ACCEPTABLE CONTROL LIMITS FOR AIR SAMPLES RESULTS > 5 x MDL | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |

Actions:

Qualify as estimated positive results (J) and nondetects (UJ) for the compound that exceeded the above criteria. For organics, only the sample and duplicate will be qualified.

If an RPD cannot be calculated because one or both of the sample results is not detected, the following actions apply:

If one sample result is not detected and the other is greater than 5x the SQL qualify (J/UJ).

If one sample value is not detected and the other is greater than 5x the SQL and the SQLs for the sample and duplicate are significantly different, use professional judgment to determine if qualification is appropriate.

If one sample value is not detected and the other is less than 5x, use professional judgment to determine if qualification is appropriate.

If both sample and duplicate results are not detected, no action is needed.

DATA REVIEW WORKSHEETS

All criteria were met X
 Criteria were not met
 and/or see below

X. INTERNAL STANDARD PERFORMANCE

The assessment of the internal standard (IS) parameter is used to assist the data reviewer in determining the condition of the analytical instrumentation.

List the internal standard area of samples which do not meet the criteria.

- * Area of +40% or -40% of the IS area in the associated calibration standard.
- * Retention time (RT) within ± 0.06 seconds of the IS area in the associated calibration standard.

| DATE | SAMPLE ID | IS OUT | IS AREA | ACCEPTABLE RANGE | ACTION |
|------|-----------|--------|---------|------------------|--------|
|------|-----------|--------|---------|------------------|--------|

Internal standard area and retention times within laboratory control limits for both samples and calibration standards except for the followings:

| | | | | | |
|----------|-----------|--------------------|--------|-----------------|-----------|
| 09/03/15 | JC2870-1A | 5-Chlorobenzene-d5 | 870406 | 150209 - 350489 | No action |
| 09/03/15 | JC2870-2 | 5-Chlorobenzene-d5 | 588889 | 150209 - 350489 | No action |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |

Actions:

1. IS actions should be applied to the compound quantitated with the out-of-control ISs

| QUALITY | IS AREA < -40% | IS AREA > + 40% |
|---------------------|----------------|-----------------|
| Positive results | J | J |
| Nondetected results | R | ACCEPT |

2. If a IS retention time varies more than 0.330 seconds, the chromatographic profile for that sample must be examined to determine if any false positive or negative exists. For shifts of a large magnitude, the reviewer may consider partial or total rejection of the data for the sample fraction.

DATA REVIEW WORKSHEETS

All criteria were met X
Criteria were not met
and/or see below

XII. SAMPLE QUANTITATION

The sample quantitation evaluation is to verify laboratory quantitation results. In the space below, please show a minimum of one sample calculation:

Cyclohexane

RF = 0.424

$$[] = (259738)(10.0)/(471861)(0.424)$$

$$= 12.98 \text{ ppbv OK}$$

DATA REVIEW WORKSHEETS

All criteria were met X
 Criteria were not met
 and/or see below

XII. QUANTITATION LIMITS

A. Dilution performed

| SAMPLE ID | DILUTION FACTOR | REASONS FOR DILUTION |
|--|-----------------|--------------------------------|
| Dilution was performed on samples by a factor of less than 3 except the following: | | |
| JC2870-1A | 59.2 | High levels of target species. |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |

B. Percent Solids

List samples which have ≤ 50 % solids

Actions:

If the % solids of a soil sample is 10-50%, estimate positive results (J) and nondetects (UJ)

If the % solids of a soil sample is < 10%, estimate positive results (J) and reject nondetects (R)

DATA REVIEW WORKSHEETS

Project Number: jc2870

Date: 08/31/2015

REVIEW OF VOLATILE ORGANIC PACKAGE

The following guidelines for evaluating volatile organics were created to delineate required validation actions. This document will assist the reviewer in using professional judgment to make more informed decision and in better serving the needs of the data users. The sample results were assessed according to USEPA data validation guidance documents in the following order of precedence: QC criteria from "Compendium Method TO-3: **METHOD FOR THE DETERMINATION OF VOLATILE ORGANIC COMPOUNDS IN AMBIENT AIR USING CRYOGENIC PRECONCENTRATION TECHNIQUES AND GAS CHROMATOGRAPHY WITH FLAME IONIZATION AND ELECTRON CAPTURE DETECTION.** The QC criteria and data validation actions listed on the data review worksheets are from the primary guidance document, unless otherwise noted.

The hardcopied (laboratory name) Accutest Laboratories data package received has been reviewed and the quality control and performance data summarized. The data review for VOCs included:

Lab. Project/SDG No.: jc2870

Sample matrix: Air

No. of Samples: 3

Trip blank No.: -

Field blank No.: -

Equipment blank No.: -

Field duplicate No.: B8-SSV1/B8-SSV1D

☒ Data Completeness

☒ Laboratory Control Spikes

☒ Holding Times

☒ Field Duplicates

☐ N/A GC/MS Tuning

☒ Calibrations

☐ N/A Internal Standard Performance

☒ Compound Identifications

☒ Blanks

☒ Compound Quantitation

☒ Surrogate Recoveries

☒ Quantitation Limits

☐ N/A Matrix Spike/Matrix Spike Duplicate

Overall Comments: Methane; Propane; n-Butane by method TO-3

Definition of Qualifiers:

J- Estimated results

U- Compound not detected

R- Rejected data

UJ- Estimated nondetect

Reviewer: Rafael Difant

Date: 09/15/2015


DATA REVIEW WORKSHEETS

DATA COMPLETENESS

MISSING INFORMATION

DATE LAB. CONTACTED

DATE RECEIVED



A large, thick, light gray diagonal line is drawn across the table area, from the top-left corner to the bottom-right corner, indicating that the data is missing or not applicable.

DATA REVIEW WORKSHEETS

All criteria were met X
Criteria were not met
and/or see below

HOLDING TIMES

The objective of this parameter is to ascertain the validity of the results based on the holding time of the sample from time of collection to the time of analysis.

Complete table for all samples and note the analysis and/or preservation not within criteria

| SAMPLE ID | DATE SAMPLED | DATE ANALYZED | pH | ACTION |
|---|--------------|---------------|----|--------|
| | | | | |
| | | | | |
| All samples analyzed within the recommended method holding time | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |

Criteria

Aqueous samples – 14 days from sample collection for preserved samples ($\text{pH} \leq 2$, 4°C), no air bubbles.

Aqueous samples – 7 days from sample collection for unpreserved samples, 4°C , no air bubbles.

Soil samples- 7 days from sample collection.

Cooler temperature (Criteria: $4 \pm 2^{\circ}\text{C}$): N/A

Actions

If the VOCs vial(s) have air bubbles, estimate positive results (J) and reject nondetects (R).

If the % solids of soil samples is 10-50%, estimate positive results (J) and nondetects (UJ)

If the % solid of soil samples is $< 10\%$, estimate positive results (J) and reject nondetects (R).

If holding times are exceeded but < 14 days beyond criteria, estimate positive results (J) and nondetects (UJ).

If holding times are exceeded but < 28 days beyond criteria, estimate positive results (J) and reject nondetects (R).

If holding times are grossly exceeded (> 28 days beyond criteria), reject all results (R).

If samples were not iced or if the ice were melted ($> 10^{\circ}\text{C}$), estimate positive results (J) and nondetects (UJ).

DATA REVIEW WORKSHEETS

All criteria were met ___N/A___
Criteria were not met see below _____

GC/MS TUNING

The assessment of the tuning results is to determine if the sample instrumentation is within the standard tuning QC limits

___N/A_ The BFB performance results were reviewed and found to be within the specified criteria.

___N/A_ BFB tuning was performed for every 24 hours of sample analysis.

If no, use professional judgment to determine whether the associated data should be accepted, qualified or rejected.

List _____ the _____ samples _____ affected:

If mass calibration is in error, all associated data are rejected.

DATA REVIEW WORKSHEETS

All criteria were met ☒
 Criteria were not met
 and/or see below _____

CALIBRATION VERIFICATION

Compliance requirements for satisfactory instrument calibration are established to ensure that the instrument is capable of producing and maintaining acceptable quantitative data.

Date of initial calibration: 09/01/15
 Dates of continuing calibration: 09/03/15; 09/04/15
 Instrument ID numbers: HP_G1530A
 Matrix/Level: Air/low

| DATE | LAB ID# | FILE | CRITERIA OUT RFs, %RSD, %D, r | COMPOUND | SAMPLES AFFECTED |
|---|------------|------|----------------------------------|----------|---------------------|
| | | | | | |
| Initial and continuing calibrations meet method specific requirements. Initial calibration retention times meet method specific requirements. | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |

Criteria

All RFs must be > 0.05 regardless of method requirements for SPCC.
 All %RSD must be $\leq 15\%$ regardless of method requirements for CCC.
 All %Ds must be $\leq 30\%$ regardless of method requirements for CCC.
 Method TO-15 does not specify criterion for the curve correlation coefficient (r). A limit for r of ≥ 0.995 has therefore been utilized as professional judgment.

Actions

If any compound has an initial RF or a continuing RF of < 0.05 , estimate positive results (J) and reject nondetects (R), regardless of method requirements.
 If any compound has a %RSD $> 15\%$, estimate positive results (J) and use professional judgment to qualify nondetects.
 If any compound has a %RSD $> 90\%$, estimate positive results (J) and reject nondetects (R).
 If any compound has a % D $> 30\%$, estimate positive results (J) and reject nondetects (R).
 If any compound has a % D $> 30\%$, estimate positive results (J) and nondetects (UJ).
 If any compound has a % D $> 90\%$, estimate positive results (J) and reject nondetects (R).
 If any compound has $r < 0.995$, estimate positive results and nondetects.

A separate worksheet should be filled for each initial curve

DATA REVIEW WORKSHEETS

All criteria were met X
Criteria were not met
and/or see below _____

V A. BLANK ANALYSIS RESULTS (Sections 1 & 2)

The assessment of the blank analysis results is to determine the existence and magnitude of contamination problems. The criteria for evaluation of blanks apply only to blanks associated with the samples, including trip, equipment, and laboratory blanks. If problems with any blanks exist, all data associated with the case must be carefully evaluated to determine whether or not there is an inherent variability in the data for the case, or if the problem is an isolated occurrence not affecting other data.

List the contamination in the blanks below. High and low levels blanks must be treated separately.

Laboratory blanks

| DATE ANALYZED | LAB ID | LEVEL/ MATRIX | COMPOUND | CONCENTRATION UNITS |
|---------------|---|------------------|----------|------------------------|
| | | | | |
| | | | | |
| | All_method_blank_meeth_method_specific_criteria | | | |
| | | | | |
| | | | | |
| | | | | |

Field/Equipment/Trip blank

[illegible]

All criteria were met X
Criteria were not met
and/or see below _____

Blank Actions

Action Levels (ALs) should be based upon the highest concentration of contaminant determined in any blank. Do not qualify any blank with another blank. The ALs for samples which have been diluted should be corrected for the sample dilution factor and/or % moisture, where applicable. No positive sample results should be reported unless the concentration of the compound in the samples exceeds the ALs:

ALs = 10x the amount of common contaminants (methylene chloride, acetone, 2-butanone, and toluene)

ALs = 5x for any other compounds.

Specific actions are as follows:

If the concentration is < sample quantitation limit (SQL) and \leq AL, report the compound as not detected (U) at the SQL.

If the concentration is \geq SQL but \leq AL, report the compound as not detected (U) at the reported concentration.

If the concentration is \geq SQL and $>$ AL, report the concentration unqualified.

Notes:

High and low level blanks must be treated separately

Compounds qualified "U" for blank contamination are still considered "hits" when qualifying for calibration criteria.

[illegible]

DATA REVIEW WORKSHEETS

All criteria were met _____
 Criteria were not met _____
 and/or see below X

SURROGATE SPIKE RECOVERIES

Laboratory performance of individual samples is established by evaluation of surrogate spike recoveries. All samples are spiked with surrogate compounds prior to sample analysis. The accuracy of the analysis is measured by the surrogate percent recovery. Since the effects of the sample matrix are frequently outside the control of the laboratory and may present relatively unique problems, the validation of data is frequently subjective and demands analytical experience and professional judgment.

List the percent recoveries (%Rs) which do not meet the criteria for surrogate recovery.

Matrix: air/solid/aqueous

| SAMPLE ID | SURROGATE COMPOUND | ACTION |
|-----------|---------------------------|----------------|
| | 1,2-DICHLOROETHANE- d4 | Toluene- d8 |
| | | 4-BFB |

 BFB used as surrogate standard; % recovery within laboratory limits except for the
 followings:

| | | |
|-----------|-------|-----------|
| JC2870-1 | 137 % | No action |
| JC2870-1D | 138 % | No action |
| | | |
| | | |
| | | |
| | | |

Note: No action taken professional judgment; % recoveries within control limits in second column.

QC Limits* (Air): LL to: 61 % UL to: 137 %

- * QC limits are laboratory in-house performance criteria, LL = lower limit, UL = upper limit.
- * If QC limits are not available, use limits of 80 – 120 % for aqueous and 70 – 130 % for solid samples.

Actions:

| QUALITY | %R < 10% | %R = 10% - LL | %R > UL |
|--------------------|----------|---------------|---------|
| Positive results | J | J | J |
| Nondetects results | R | UJ | Accept |

Surrogate action should be applied:

If one or more surrogate in the VOC fraction is out of specification, but has a recovery of > 10%.

If any one surrogate in a fraction shows < 10 % recovery.

All criteria were met _____
 Criteria were not met _____
 and/or see below ___N/A___

VII. A MATRIX SPIKE/MATRIX SPIKE DUPLICATE (MS/MSD)

This data is generated to determine long term precision and accuracy in the analytical method for various matrices. This data alone cannot be used to evaluate the precision and accuracy of individual samples. If any % R in the MS or MSD falls outside the designated range, the reviewer should determine if there are matrix effects, i.e. LCS data are within the QC limits but MS/MSD data are outside QC limit.

1. MS/MSD Recoveries and Precision Criteria

The laboratory should use one MS and a duplicate analysis of an unspiked field sample if target analytes are expected in the sample. If target analytes are not expected, MS/MSD should be analyzed.

List the %Rs, RPD of the compounds which do not meet the criteria.

Sample ID: _____ Matrix/Level: _____

| MS OR MSD | COMPOUND | % R | RPD | QC LIMITS | ACTION |
|--|----------|-----|-----|-----------|--------|
| _____MS/MSD are not required as part of the method; blank spike used to assess accuracy_____ | | | | | |
| _____ | | | | | |

- * QC limits are laboratory in-house performance criteria, LL = lower limit, UL = upper limit.
- * If QC limits are not available, use limits of 70 – 130 %.

Actions:

| QUALITY | %R < LL | %R > UL |
|--------------------|---------|---------|
| Positive results | J | J |
| Nondetects results | R | Accept |

MS/MSD criteria apply only to the unspiked sample, its dilutions, and the associated MS/MSD samples:

If the % R for the affected compounds were < LL (or 70 %), qualify positive results (J) and nondetects (JJ).

If the % R for the affected compounds were > UL (or 130 %), only qualify positive results (J).

If 25 % or more of all MS/MSD %R were < LL (or 70 %) or if two or more MS/MSD %Rs were < 10%, qualify all positive results (J) and reject nondetects (R).

A separate worksheet should be used for each MS/MSD pair.

All criteria were met _____
Criteria were not met _____
and/or see below ___N/A___

MS/MSD – Unspiked Compounds

If all target analytes were spiked in the MS/MSD, this review element is not applicable.

Sample ID: _____ Matrix/Level/Unit: _____

[illegible]

* If the % RSD > 50, qualify the positive result in the unspiked samples as estimated (J).
* If the % RSD is not calculated (NC) due to nondetected value, use professional judgment to qualify the data.

DATA REVIEW WORKSHEETS

All criteria were met X
 Criteria were not met
 and/or see below

VIII. LABORATORY CONTROL SAMPLE (LCS) ANALYSIS

This data is generated to determine accuracy of the analytical method for various matrices.

1. LCS Recoveries Criteria

Where LCS spiked with the same analyte at the same concentrations as the MS/MSD?
 Yes or No. If no make note in data review memo.

List the %R of compounds which do not meet the criteria

| LCS ID | COMPOUND | % R | QC LIMIT |
|---|----------|-----|----------|
| LCS (Blank spike) analyzed in this data package, recoveries and RPD within laboratory control limits. | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |
| | | | |

- * QC limits are laboratory in-house performance criteria, LL = lower limit, UL = upper limit.
- * If QC limits are not available, use limits of 70 – 130 %.

Actions:

| QUALITY | %R < LL | %R > UL |
|--------------------|---------|---------|
| Positive results | J | J |
| Nondetects results | R | Accept |

All analytes in the associated sample results are qualified for the following criteria.

If 25 % of the LCS recoveries were < LL (or 70 %), qualify all positive results (j) and reject nondetects (R).

If two or more LCS were below 10 %, qualify all positive results as (J) and reject nondetects (R).

2. Frequency Criteria:

Where LCS analyzed at the required frequency and for each matrix? Yes or No.

If no, the data may be affected. Use professional judgment to determine the severity of the effect and qualify data accordingly. Discuss any actions below and list the samples affected.

DATA REVIEW WORKSHEETS

All criteria were met X
 Criteria were not met
 and/or see below

IX. LABORATORY DUPLICATE PRECISION

Sample IDs: JC2870-1/-1_DUP
 Sample IDs: JC3132-1/-1_DUP

Matrix: Air
 Matrix: Air

Field duplicate samples may be taken and analyzed as an indication of overall precision. These analyses measure both field and lab precision; therefore, the results may have more variability than laboratory duplicates which only laboratory performance. It is also expected that soil duplicate results will have a greater variance than water matrices due to difficulties associated with collecting identical field duplicate samples.

The project QAPP should be reviewed for project-specific information.

Suggested criteria: RPD \pm 25% for air samples. If both samples and duplicate are <5 SQL, the RPD criteria is doubled.

| COMPOUND | SQL | SAMPLE CONC. | DUPLICATE CONC. | RPD | ACTION |
|---|-----|--------------|-----------------|-----|--------|
| | | | | | |
| | | | | | |
| RPD within the method performance criteria. | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |

Actions:

Qualify as estimated positive results (J) and nondetects (UJ) for the compound that exceeded the above criteria. For organics, only the sample and duplicate will be qualified.

If an RPD cannot be calculated because one or both of the sample results is not detected, the following actions apply:

If one sample result is not detected and the other is greater than 5x the SQL qualify (J/UJ).

If one sample value is not detected and the other is greater than 5x the SQL and the SQLs for the sample and duplicate are significantly different, use professional judgment to determine if qualification is appropriate.

If one sample value is not detected and the other is less than 5x, use professional judgment to determine if qualification is appropriate.

If both sample and duplicate results are not detected, no action is needed.

DATA REVIEW WORKSHEETS

All criteria were met X
 Criteria were not met
 and/or see below _____

IX. FIELD DUPLICATE PRECISION

Sample IDs: JC2870-2/JC2870-3

Matrix: Air

Field duplicate samples may be taken and analyzed as an indication of overall precision. These analyses measure both field and lab precision; therefore, the results may have more variability than laboratory duplicates which only laboratory performance. It is also expected that soil duplicate results will have a greater variance than water matrices due to difficulties associated with collecting identical field duplicate samples.

The project QAPP should be reviewed for project-specific information.

Suggested criteria: RPD \pm 50% for air samples. If both samples and duplicate are <5 SQL, the RPD criteria is doubled.

| COMPOUND | SQL | SAMPLE CONC. | DUPLICATE CONC. | RPD | ACTION |
|----------|-----------|--------------|-----------------|-----|--|
| Methane | 0.33/0.51 | 121 | 197 | 48 | Within generally acceptable control limits |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |

Actions:

Qualify as estimated positive results (J) and nondetects (UJ) for the compound that exceeded the above criteria. For organics, only the sample and duplicate will be qualified.

If an RPD cannot be calculated because one or both of the sample results is not detected, the following actions apply:

If one sample result is not detected and the other is greater than 5x the SQL qualify (J/UJ).

If one sample value is not detected and the other is greater than 5x the SQL and the SQLs for the sample and duplicate are significantly different, use professional judgment to determine if qualification is appropriate.

If one sample value is not detected and the other is less than 5x, use professional judgment to determine if qualification is appropriate.

If both sample and duplicate results are not detected, no action is needed.

DATA REVIEW WORKSHEETS

All criteria were met __N/A__
Criteria were not met
and/or see below ____

X. INTERNAL STANDARD PERFORMANCE

The assessment of the internal standard (IS) parameter is used to assist the data reviewer in determining the condition of the analytical instrumentation.

List the internal standard area of samples which do not meet the criteria.

- * Area of +40% or -40% of the IS area in the associated calibration standard.
- * Retention time (RT) within ± 0.06 seconds of the IS area in the associated calibration standard.

| DATE | SAMPLE ID | IS OUT | IS AREA | ACCEPTABLE RANGE | ACTION |
|------|-----------|--------|---------|------------------|--------|
|------|-----------|--------|---------|------------------|--------|

Internal standard not used for calibration; samples quantitated by the external standard method.

| |
|--|
| |
| |
| |
| |
| |

Actions:

1. IS actions should be applied to the compound quantitated with the out-of-control ISs

| QUALITY | IS AREA < -40% | | IS AREA > + 40% |
|---------------------|----------------|--|-----------------|
| Positive results | J | | J |
| Nondetected results | R | | ACCEPT |

2. If a IS retention time varies more than 0.330 seconds, the chromatographic profile for that sample must be examined to determine if any false positive or negative exists. For shifts of a large magnitude, the reviewer may consider partial or total rejection of the data for the sample fraction.

DATA REVIEW WORKSHEETS

All criteria were met X
Criteria were not met
and/or see below

XII. SAMPLE QUANTITATION

The sample quantitation evaluation is to verify laboratory quantitation results. In the space below, please show a minimum of one sample calculation:

JC2870-1

Methane RF = 583900

$$[] = (75908746108)/583900$$

$$= 1.3 \times 10^5 \text{ ppmv OK}$$

A. Dilution performed

[illegible]

List samples which have $\leq 50\%$ solids

A thick, dark gray diagonal line with a fine, cross-hatched texture. It runs from the bottom-left towards the top-right, crossing the horizontal lines of the page.

If the % solids of a soil sample is $< 10\%$, estimate positive results (J) and reject nondetects (R)