

CERTIFICATION

SDG No:	MC46423	Laboratory:	Accutest, Massachusetts
Site:	BMS, Building 5 Area, PR Humacao, PR	Matrix:	Groundwater

SUMMARY: Groundwater samples (Table 1) were collected on the BMSMC facility – Building 5 Area. The BMSMC facility is located in Humacao, PR. Samples were taken July 13-16, 2016 and were analyzed in Accutest Laboratory of Marlborough, Massachusetts that reported the data under SDG No.: MC46423. Results were validated using the following quality control criteria of the methods employed (MADEP VPH and MAPED EPH, Massachusetts Department of Environmental Protection, 2004) and the latest validation guidelines (July, 2015) of the EPA Hazardous Waste Support Section. The analyses performed are shown in Table 1. Individual data review worksheets are enclosed for each target analyte group. The data sample organic data samples summary form shows for analytes results that were qualified.

In summary the results are valid and can be used for decision taking purposes.

Table 1. Samples analyzed and analysis performed

SAMPLE ID	SAMPLE DESCRIPTION	MATRIX	ANALYSIS PERFORMED
MC46423-1	S-33	Groundwater	Volatiles TPHC Ranges
MC46423-1A	S-33	Groundwater	Extractable TPHC Ranges
MC46423-2	S-34	Groundwater	Volatiles TPHC Ranges
MC46423-2A	S-34	Groundwater	Extractable TPHC Ranges
MC46423-3	G-1R3	Groundwater	Volatiles TPHC Ranges
MC46423-3A	G-1R3	Groundwater	Extractable TPHC Ranges
MC46423-4	E-1R	Groundwater	Volatiles TPHC Ranges
MC46423-4A	E-1R	Groundwater	Extractable TPHC Ranges
MC46423-5	D-1R	Groundwater	Volatiles TPHC Ranges
MC46423-5A	D-1R	Groundwater	Extractable TPHC Ranges
MC46423-6	MW-19	Groundwater	Volatile TPHC Ranges
MC46423-6A	MW-19	Groundwater	Extractable TPHC Ranges

SAMPLE ID	SAMPLE DESCRIPTION	MATRIX	ANALYSIS PERFORMED
MC46423-7	MW-22S	Groundwater	Volatile TPHC Ranges
MC46423-7A	MW-22S	Groundwater	Extractable TPHC Ranges

Reviewer Name: Rafael Infante
Chemist License 1888

Signature:

Rafael Infante

Date:

July 16, 2016



SGS Accutest

Report of Analysis

Page 1 of 1

Client Sample ID:	S-33	Date Sampled:	06/13/16
Lab Sample ID:	MC46423-1	Date Received:	06/17/16
Matrix:	AQ - Ground Water	Percent Solids:	n/a
Method:	MADEP VPH REV 1.1		
Project:	BMSMC, Building 5 Area, Puerto Rico		

Run #	File ID	DF	Analyzed	By	Prep Date	Prep Batch	Analytical Batch
Run #1	WX77172.D	1	06/20/16	AF	n/a	n/a	GWX3796
Run #2							

Run #	Purge Volume
Run #1	5.0 ml
Run #2	

Volatile TPHC Ranges

CAS No.	Compound	Result	RL	MDL	Units	Q
	C5- C8 Aliphatics (Unadj.)	30.9	50	25	ug/l	J
	C9- C12 Aliphatics (Unadj.)	72.3	50	25	ug/l	
	C9- C10 Aromatics (Unadj.)	50.0	50	25	ug/l	
	C5- C8 Aliphatics	ND	50	25	ug/l	
	C9- C12 Aliphatics	ND	50	25	ug/l	

CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limits
	2,3,4-Trifluorotoluene	98%		70-130%
	2,3,4-Trifluorotoluene	101%		70-130%



ND = Not detected MDL = Method Detection Limit
 RL = Reporting Limit
 E = Indicates value exceeds calibration range

J = Indicates an estimated value
 B = Indicates analyte found in associated method blank
 N = Indicates presumptive evidence of a compound

SGS Accutest

Report of Analysis

Page 1 of 1

Client Sample ID:	S-33	Date Sampled:	06/13/16
Lab Sample ID:	MC46423-1A	Date Received:	06/17/16
Matrix:	AQ - Ground Water	Percent Solids:	n/a
Method:	MADEP EPH REV 1.1 SW846 3510C		
Project:	BMSMC, Building 5 Area, Puerto Rico		

Run #	File ID	DF	Analyzed	By	Prep Date	Prep Batch	Analytical Batch
Run #1	DE14765.D	1	06/30/16	TA	06/27/16	OP47988	GDE820
Run #2							

Run #	Initial Volume	Final Volume
Run #1	950 ml	2.0 ml
Run #2		

Extractable TPHC Ranges

CAS No.	Compound	Result	RL	MDL	Units	Q
	C11-C22 Aromatics (Unadj.)	33.2	110	30	ug/l	JB
	C9-C18 Aliphatics	20.6	110	18	ug/l	J
	C19-C36 Aliphatics	32.9	110	29	ug/l	J
	C11-C22 Aromatics	32.6	110	30	ug/l	JB

CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limits
84-15-1	o-Terphenyl	45%		40-140%
321-60-8	2-Fluorobiphenyl	70%		40-140%
3386-33-2	1-Chlorooctadecane	43%		40-140%
580-13-2	2-Bromonaphthalene	77%		40-140%



ND = Not detected MDL = Method Detection Limit
 RL = Reporting Limit
 E = Indicates value exceeds calibration range

J = Indicates an estimated value
 B = Indicates analyte found in associated method blank
 N = Indicates presumptive evidence of a compound

SGS Accutest

Report of Analysis

Page 1 of 1

Client Sample ID:	S-34	Date Sampled:	06/13/16
Lab Sample ID:	MC46423-2	Date Received:	06/17/16
Matrix:	AQ - Ground Water	Percent Solids:	n/a
Method:	MADEP VPH REV 1.1		
Project:	BMSMC, Building 5 Area, Puerto Rico		

Run #	File ID	DF	Analyzed	By	Prep Date	Prep Batch	Analytical Batch
Run #1	WX77173.D	1	06/20/16	AF	n/a	n/a	GWX3796
Run #2							

Run #	Purge Volume
Run #1	5.0 ml
Run #2	

Volatile TPHC Ranges

CAS No.	Compound	Result	RL	MDL	Units	Q
	C5- C8 Aliphatics (Unadj.)	ND	50	25	ug/l	
	C9- C12 Aliphatics (Unadj.)	ND	50	25	ug/l	
	C9- C10 Aromatics (Unadj.)	ND	50	25	ug/l	
	C5- C8 Aliphatics	ND	50	25	ug/l	
	C9- C12 Aliphatics	ND	50	25	ug/l	

CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limits
	2,3,4-Trifluorotoluene	94%		70-130%
	2,3,4-Trifluorotoluene	96%		70-130%



ND = Not detected MDL = Method Detection Limit
 RL = Reporting Limit
 E = Indicates value exceeds calibration range

J = Indicates an estimated value
 B = Indicates analyte found in associated method blank
 N = Indicates presumptive evidence of a compound

SGS Accutest

Report of Analysis

Page 1 of 1

Client Sample ID:	S-34	Date Sampled:	06/13/16
Lab Sample ID:	MC46423-2A	Date Received:	06/17/16
Matrix:	AQ - Ground Water	Percent Solids:	n/a
Method:	MADEP EPH REV 1.1 SW846 3510C		
Project:	BMSMC, Building 5 Area, Puerto Rico		

Run #	File ID	DF	Analyzed	By	Prep Date	Prep Batch	Analytical Batch
Run #1	DE14766.D	1	06/30/16	TA	06/27/16	OP47988	GDE820
Run #2							

Run #	Initial Volume	Final Volume
Run #1	920 ml	2.0 ml
Run #2		

Extractable TPHC Ranges

CAS No.	Compound	Result	RL	MDL	Units	Q
	C11-C22 Aromatics (Unadj.)	34.2	110	31	ug/l	JB
	C9-C18 Aliphatics	22.7	110	18	ug/l	J
	C19-C36 Aliphatics	40.9	110	29	ug/l	J
	C11-C22 Aromatics	34.2	110	31	ug/l	JB

CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limits
84-15-1	o-Terphenyl	63%		40-140%
321-60-8	2-Fluorobiphenyl	71%		40-140%
3386-33-2	1-Chlorooctadecane	63%		40-140%
580-13-2	2-Bromonaphthalene	78%		40-140%



ND = Not detected MDL = Method Detection Limit
 RL = Reporting Limit
 E = Indicates value exceeds calibration range

J = Indicates an estimated value
 B = Indicates analyte found in associated method blank
 N = Indicates presumptive evidence of a compound

SGS Accutest

Report of Analysis

Page 1 of 1

Client Sample ID:	G-1R3	Date Sampled:	06/15/16
Lab Sample ID:	MC46423-3	Date Received:	06/17/16
Matrix:	AQ - Ground Water	Percent Solids:	n/a
Method:	MADEP VPH REV 1.1		
Project:	BMSMC, Building 5 Area, Puerto Rico		

Run #	File ID	DF	Analyzed	By	Prep Date	Prep Batch	Analytical Batch
Run #1	WX77174.D	1	06/20/16	AF	n/a	n/a	GWX3796
Run #2	WX77180.D	100	06/20/16	AF	n/a	n/a	GWX3796

Run #	Purge Volume
Run #1	5.0 ml
Run #2	5.0 ml

Volatile TPHC Ranges

CAS No.	Compound	Result	RL	MDL	Units	Q
	C5- C8 Aliphatics (Unadj.)	167	50	25	ug/l	
	C9- C12 Aliphatics (Unadj.)	63100 ^a	5000	2500	ug/l	
	C9- C10 Aromatics (Unadj.)	112	50	25	ug/l	
	C5- C8 Aliphatics	67.1	50	25	ug/l	
	C9- C12 Aliphatics	753	50	25	ug/l	
CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limits		
	2,3,4-Trifluorotoluene	106%	90%	70-130%		
	2,3,4-Trifluorotoluene	108%	94%	70-130%		

(a) Result is from Run# 2



ND = Not detected MDL = Method Detection Limit
 RL = Reporting Limit
 E = Indicates value exceeds calibration range

J = Indicates an estimated value
 B = Indicates analyte found in associated method blank
 N = Indicates presumptive evidence of a compound

SGS Accutest

Report of Analysis

Page 1 of 1

Client Sample ID:	G-1R3	Date Sampled:	06/15/16
Lab Sample ID:	MC46423-3A	Date Received:	06/17/16
Matrix:	AQ - Ground Water	Percent Solids:	n/a
Method:	MADEP EPH REV 1.1 SW846 3510C		
Project:	BMSMC, Building 5 Area, Puerto Rico		

Run #	File ID	DF	Analyzed	By	Prep Date	Prep Batch	Analytical Batch
Run #1	DE14767.D	1	06/30/16	TA	06/27/16	OP47988	GDE820
Run #2							

Run #	Initial Volume	Final Volume
Run #1	935 ml	2.0 ml
Run #2		

Extractable TPHC Ranges

CAS No.	Compound	Result	RL	MDL	Units	Q
	C11-C22 Aromatics (Unadj.)	36.4	110	31	ug/l	JB
	C9-C18 Aliphatics	25.8	110	18	ug/l	J
	C19-C36 Aliphatics	96.7	110	29	ug/l	J
	C11-C22 Aromatics	36.4	110	31	ug/l	JB

CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limits
84-15-1	o-Terphenyl	51%		40-140%
321-60-8	2-Fluorobiphenyl	72%		40-140%
3386-33-2	1-Chlorooctadecane	60%		40-140%
580-13-2	2-Bromonaphthalene	78%		40-140%



ND = Not detected MDL = Method Detection Limit
 RL = Reporting Limit
 E = Indicates value exceeds calibration range

J = Indicates an estimated value
 B = Indicates analyte found in associated method blank
 N = Indicates presumptive evidence of a compound

SGS Accutest

Report of Analysis

Page 1 of 1

Client Sample ID: E-1R
 Lab Sample ID: MC46423-4
 Matrix: AQ - Ground Water
 Method: MADEP VPH REV 1.1
 Project: BMSMC, Building 5 Area, Puerto Rico

Date Sampled: 06/15/16
 Date Received: 06/17/16
 Percent Solids: n/a

Run #	File ID	DF	Analyzed	By	Prep Date	Prep Batch	Analytical Batch
Run #1	WX77179.D	1	06/20/16	AF	n/a	n/a	GWX3796
Run #2							

Run #	Purge Volume
Run #1	5.0 ml
Run #2	

Volatile TPHC Ranges

CAS No.	Compound	Result	RL	MDL	Units	Q
	C5- C8 Aliphatics (Unadj.)	36.4	50	25	ug/l	J
	C9- C12 Aliphatics (Unadj.)	55.4	50	25	ug/l	
	C9- C10 Aromatics (Unadj.)	ND	50	25	ug/l	
	C5- C8 Aliphatics	27.7	50	25	ug/l	J
	C9- C12 Aliphatics	ND	50	25	ug/l	

CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limits
	2,3,4-Trifluorotoluene	91%		70-130%
	2,3,4-Trifluorotoluene	94%		70-130%



ND = Not detected MDL = Method Detection Limit
 RL = Reporting Limit
 E = Indicates value exceeds calibration range

J = Indicates an estimated value
 B = Indicates analyte found in associated method blank
 N = Indicates presumptive evidence of a compound

SGS Accutest

Report of Analysis

Page 1 of 1

Client Sample ID:	E-1R	Date Sampled:	06/15/16
Lab Sample ID:	MC46423-4A	Date Received:	06/17/16
Matrix:	AQ - Ground Water	Percent Solids:	n/a
Method:	MADEP EPH REV 1.1 SW846 3510C		
Project:	BMSMC, Building 5 Area, Puerto Rico		

Run #	File ID	DF	Analyzed	By	Prep Date	Prep Batch	Analytical Batch
Run #1	DE14768.D	1	06/30/16	TA	06/27/16	OP47988	GDE820
Run #2							

Run #	Initial Volume	Final Volume
Run #1	970 ml	2.0 ml
Run #2		

Extractable TPHC Ranges

CAS No.	Compound	Result	RL	MDL	Units	Q
	C11-C22 Aromatics (Unadj.)	60.6	100	30	ug/l	JB
	C9-C18 Aliphatics	22.7	100	17	ug/l	J
	C19-C36 Aliphatics	46.6	100	28	ug/l	J
	C11-C22 Aromatics	32.4	100	30	ug/l	JB

CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limits
84-15-1	o-Terphenyl	61%		40-140%
321-60-8	2-Fluorobiphenyl	67%		40-140%
3386-33-2	1-Chlorooctadecane	64%		40-140%
580-13-2	2-Bromonaphthalene	74%		40-140%



ND = Not detected MDL = Method Detection Limit
 RL = Reporting Limit
 E = Indicates value exceeds calibration range

J = Indicates an estimated value
 B = Indicates analyte found in associated method blank
 N = Indicates presumptive evidence of a compound

SGS Accutest

Report of Analysis

Page 1 of 1

Client Sample ID:	D-1R	Date Sampled:	06/15/16
Lab Sample ID:	MC46423-5	Date Received:	06/17/16
Matrix:	AQ - Ground Water	Percent Solids:	n/a
Method:	MADEP VPH REV 1.1		
Project:	BMSMC, Building 5 Area, Puerto Rico		

Run #	File ID	DF	Analyzed	By	Prep Date	Prep Batch	Analytical Batch
Run #1	WX77176.D	1	06/20/16	AF	n/a	n/a	GWX3796
Run #2							

Run #	Purge Volume
Run #1	5.0 ml
Run #2	

Volatile TPHC Ranges

CAS No.	Compound	Result	RL	MDL	Units	Q
	C5- C8 Aliphatics (Unadj.)	ND	50	25	ug/l	
	C9- C12 Aliphatics (Unadj.)	ND	50	25	ug/l	
	C9- C10 Aromatics (Unadj.)	ND	50	25	ug/l	
	C5- C8 Aliphatics	ND	50	25	ug/l	
	C9- C12 Aliphatics	ND	50	25	ug/l	
CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limits		
	2,3,4-Trifluorotoluene	90%		70-130%		
	2,3,4-Trifluorotoluene	93%		70-130%		



ND = Not detected MDL = Method Detection Limit
 RL = Reporting Limit
 E = Indicates value exceeds calibration range

J = Indicates an estimated value
 B = Indicates analyte found in associated method blank
 N = Indicates presumptive evidence of a compound

SGS Accutest

Report of Analysis

Page 1 of 1

Client Sample ID:	D-1R	Date Sampled:	06/15/16
Lab Sample ID:	MC46423-5A	Date Received:	06/17/16
Matrix:	AQ - Ground Water	Percent Solids:	n/a
Method:	MADEP EPH REV 1.1 SW846 3510C		
Project:	BMSMC, Building 5 Area, Puerto Rico		

Run #	File ID	DF	Analyzed	By	Prep Date	Prep Batch	Analytical Batch
Run #1	DE14769.D	1	06/30/16	TA	06/27/16	OP47988	GDE820
Run #2							

Run #	Initial Volume	Final Volume
Run #1	940 ml	2.0 ml
Run #2		

Extractable TPHC Ranges

CAS No.	Compound	Result	RL	MDL	Units	Q
	C11-C22 Aromatics (Unadj.)	ND	110	30	ug/l	
	C9-C18 Aliphatics	20.8	110	18	ug/l	J
	C19-C36 Aliphatics	33.7	110	29	ug/l	J
	C11-C22 Aromatics	ND	110	30	ug/l	

CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limits
84-15-1	o-Terphenyl	42%		40-140%
321-60-8	2-Fluorobiphenyl	67%		40-140%
3386-33-2	1-Chlorooctadecane	43%		40-140%
580-13-2	2-Bromonaphthalene	74%		40-140%



ND = Not detected MDL = Method Detection Limit
 RL = Reporting Limit
 E = Indicates value exceeds calibration range

J = Indicates an estimated value
 B = Indicates analyte found in associated method blank
 N = Indicates presumptive evidence of a compound

SGS Accutest

Report of Analysis

Page 1 of 1

Client Sample ID:	MW-19	Date Sampled:	06/16/16
Lab Sample ID:	MC46423-6	Date Received:	06/17/16
Matrix:	AQ - Ground Water	Percent Solids:	n/a
Method:	MADEP VPH REV 1.1		
Project:	BMSMC, Building 5 Area, Puerto Rico		

Run #	File ID	DF	Analyzed	By	Prep Date	Prep Batch	Analytical Batch
Run #1	WX77177.D	1	06/20/16	AF	n/a	n/a	GWX3796
Run #2	WX77181.D	50	06/20/16	AF	n/a	n/a	GWX3796

Run #	Purge Volume
Run #1	5.0 ml
Run #2	5.0 ml

Volatile TPHC Ranges

CAS No.	Compound	Result	RL	MDL	Units	Q
	C5- C8 Aliphatics (Unadj.)	27.2	50	25	ug/l	J
	C9- C12 Aliphatics (Unadj.)	19500 ^a	2500	1300	ug/l	
	C9- C10 Aromatics (Unadj.)	100	50	25	ug/l	
	C5- C8 Aliphatics	26.3	50	25	ug/l	J
	C9- C12 Aliphatics	65.3	50	25	ug/l	
CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limits		
	2,3,4-Trifluorotoluene	93%	91%	70-130%		
	2,3,4-Trifluorotoluene	96%	95%	70-130%		

(a) Result is from Run# 2



ND = Not detected MDL = Method Detection Limit
 RL = Reporting Limit
 E = Indicates value exceeds calibration range

J = Indicates an estimated value
 B = Indicates analyte found in associated method blank
 N = Indicates presumptive evidence of a compound

SGS Accutest

Report of Analysis

Page 1 of 1

Client Sample ID:	MW-19	Date Sampled:	06/16/16
Lab Sample ID:	MC46423-6A	Date Received:	06/17/16
Matrix:	AQ - Ground Water	Percent Solids:	n/a
Method:	MADEP EPH REV 1.1 SW846 3510C		
Project:	BMSMC, Building 5 Area, Puerto Rico		

Run #	File ID	DF	Analyzed	By	Prep Date	Prep Batch	Analytical Batch
Run #1	DE14770.D	1	06/30/16	TA	06/27/16	OP47988	GDE820
Run #2							

Run #	Initial Volume	Final Volume
Run #1	940 ml	2.0 ml
Run #2		

Extractable TPHC Ranges

CAS No.	Compound	Result	RL	MDL	Units	Q
	C11-C22 Aromatics (Unadj.)	85.3	110	30	ug/l	JB
	C9-C18 Aliphatics	72.7	110	18	ug/l	J
	C19-C36 Aliphatics	37.3	110	29	ug/l	J
	C11-C22 Aromatics	78.4	110	30	ug/l	JB

CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limits
84-15-1	o-Terphenyl	51%		40-140%
321-60-8	2-Fluorobiphenyl	64%		40-140%
3386-33-2	1-Chlorooctadecane	48%		40-140%
580-13-2	2-Bromonaphthalene	67%		40-140%



ND = Not detected MDL = Method Detection Limit
 RL = Reporting Limit
 E = Indicates value exceeds calibration range

J = Indicates an estimated value
 B = Indicates analyte found in associated method blank
 N = Indicates presumptive evidence of a compound

SGS Accutest

Report of Analysis

Page 1 of 1

Client Sample ID:	MW-22S	Date Sampled:	06/16/16
Lab Sample ID:	MC46423-7	Date Received:	06/17/16
Matrix:	AQ - Ground Water	Percent Solids:	n/a
Method:	MADEP VPH REV 1.1		
Project:	BMSMC, Building 5 Area, Puerto Rico		

Run #	File ID	DF	Analyzed	By	Prep Date	Prep Batch	Analytical Batch
Run #1	WX77182.D	1	06/20/16	AF	n/a	n/a	GWX3796
Run #2							

Run #	Purge Volume
Run #1	5.0 ml
Run #2	

Volatile TPHC Ranges

CAS No.	Compound	Result	RL	MDL	Units	Q
	C5- C8 Aliphatics (Unadj.)	ND	50	25	ug/l	
	C9- C12 Aliphatics (Unadj.)	ND	50	25	ug/l	
	C9- C10 Aromatics (Unadj.)	ND	50	25	ug/l	
	C5- C8 Aliphatics	ND	50	25	ug/l	
	C9- C12 Aliphatics	ND	50	25	ug/l	
CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limits		
	2,3,4-Trifluorotoluene	90%		70-130%		
	2,3,4-Trifluorotoluene	94%		70-130%		



ND = Not detected MDL = Method Detection Limit
 RL = Reporting Limit
 E = Indicates value exceeds calibration range

J = Indicates an estimated value
 B = Indicates analyte found in associated method blank
 N = Indicates presumptive evidence of a compound

SGS Accutest

Report of Analysis

Page 1 of 1

Client Sample ID:	MW-22S	Date Sampled:	06/16/16
Lab Sample ID:	MC46423-7A	Date Received:	06/17/16
Matrix:	AQ - Ground Water	Percent Solids:	n/a
Method:	MADEP EPH REV 1.1 SW846 3510C		
Project:	BMSMC, Building 5 Area, Puerto Rico		

	File ID	DF	Analyzed	By	Prep Date	Prep Batch	Analytical Batch
Run #1	DE14771.D	1	06/30/16	TA	06/27/16	OP47988	GDE820
Run #2							

	Initial Volume	Final Volume
Run #1	920 ml	2.0 ml
Run #2		

Extractable TPHC Ranges

CAS No.	Compound	Result	RL	MDL	Units	Q
	C11-C22 Aromatics (Unadj.)	39.1	110	31	ug/l	JB
	C9-C18 Aliphatics	19.9	110	18	ug/l	J
	C19-C36 Aliphatics	59.6	110	29	ug/l	J
	C11-C22 Aromatics	39.1	110	31	ug/l	JB

CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limits
84-15-1	o-Terphenyl	60%		40-140%
321-60-8	2-Fluorobiphenyl	67%		40-140%
3386-33-2	1-Chlorooctadecane	65%		40-140%
580-13-2	2-Bromonaphthalene	74%		40-140%



ND = Not detected MDL = Method Detection Limit
 RL = Reporting Limit
 E = Indicates value exceeds calibration range

J = Indicates an estimated value
 B = Indicates analyte found in associated method blank
 N = Indicates presumptive evidence of a compound

SGS

ACCUTEST

CHAIN OF CUSTODY

SGS Accutest - Dept. New England
2334 River 190 - Daytonville, MA 01828
TEL 732-329-0200 FAX 732-329-3499/1450
www.accutest.com Marlborough MA

809960562870

PAGE 1 OF 1

Batch Order Control # MC46423
SGS Account # 6017 from NJ

Client / Reporting Information		Project Information		Requested Analysis (see TEST CODE sheet)												Matrix Codes																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																		
Company Name Anderson Mulholland Assoc. Inc.		Project Name BMS																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																
Address 2700 Westchester Purchase NY		City Humacao PR																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																
Project Contact Larry Taylor		Project #																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																
Phone # 914-251-0400		Client Purchase Order #																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																
Signature (Name) Nestor M. Rivera		Project Manager																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																
Field ID / Point of Collection		Collection		Number of (pretreated) bottles												LAB USE ONLY																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																		
		Date Time		Sampled By		Matrix		# of bottles		PO	HAH1	HAH2	HAH3	HAH4	HAH5	HAH6	HAH7	HAH8	HAH9	HAH10	HAH11	HAH12	HAH13	HAH14	HAH15	HAH16	HAH17	HAH18	HAH19	HAH20	HAH21	HAH22	HAH23	HAH24	HAH25	HAH26	HAH27	HAH28	HAH29	HAH30	HAH31	HAH32	HAH33	HAH34	HAH35	HAH36	HAH37	HAH38	HAH39	HAH40	HAH41	HAH42	HAH43	HAH44	HAH45	HAH46	HAH47	HAH48	HAH49	HAH50	HAH51	HAH52	HAH53	HAH54	HAH55	HAH56	HAH57	HAH58	HAH59	HAH60	HAH61	HAH62	HAH63	HAH64	HAH65	HAH66	HAH67	HAH68	HAH69	HAH70	HAH71	HAH72	HAH73	HAH74	HAH75	HAH76	HAH77	HAH78	HAH79	HAH80	HAH81	HAH82	HAH83	HAH84	HAH85	HAH86	HAH87	HAH88	HAH89	HAH90	HAH91	HAH92	HAH93	HAH94	HAH95	HAH96	HAH97	HAH98	HAH99	HAH100	HAH101	HAH102	HAH103	HAH104	HAH105	HAH106	HAH107	HAH108	HAH109	HAH110	HAH111	HAH112	HAH113	HAH114	HAH115	HAH116	HAH117	HAH118	HAH119	HAH120	HAH121	HAH122	HAH123	HAH124	HAH125	HAH126	HAH127	HAH128	HAH129	HAH130	HAH131	HAH132	HAH133	HAH134	HAH135	HAH136	HAH137	HAH138	HAH139	HAH140	HAH141	HAH142	HAH143	HAH144	HAH145	HAH146	HAH147	HAH148	HAH149	HAH150	HAH151	HAH152	HAH153	HAH154	HAH155	HAH156	HAH157	HAH158	HAH159	HAH160	HAH161	HAH162	HAH163	HAH164	HAH165	HAH166	HAH167	HAH168	HAH169	HAH170	HAH171	HAH172	HAH173	HAH174	HAH175	HAH176	HAH177	HAH178	HAH179	HAH180	HAH181	HAH182	HAH183	HAH184	HAH185	HAH186	HAH187	HAH188	HAH189	HAH190	HAH191	HAH192	HAH193	HAH194	HAH195	HAH196	HAH197	HAH198	HAH199	HAH200	HAH201	HAH202	HAH203	HAH204	HAH205	HAH206	HAH207	HAH208	HAH209	HAH210	HAH211	HAH212	HAH213	HAH214	HAH215	HAH216	HAH217	HAH218	HAH219	HAH220	HAH221	HAH222	HAH223	HAH224	HAH225	HAH226	HAH227	HAH228	HAH229	HAH230	HAH231	HAH232	HAH233	HAH234	HAH235	HAH236	HAH237	HAH238	HAH239	HAH240	HAH241	HAH242	HAH243	HAH244	HAH245	HAH246	HAH247	HAH248	HAH249	HAH250	HAH251	HAH252	HAH253	HAH254	HAH255	HAH256	HAH257	HAH258	HAH259	HAH260	HAH261	HAH262	HAH263	HAH264	HAH265	HAH266	HAH267	HAH268	HAH269	HAH270	HAH271	HAH272	HAH273	HAH274	HAH275	HAH276	HAH277	HAH278	HAH279	HAH280	HAH281	HAH282	HAH283	HAH284	HAH285	HAH286	HAH287	HAH288	HAH289	HAH290	HAH291	HAH292	HAH293	HAH294	HAH295	HAH296	HAH297	HAH298	HAH299	HAH300	HAH301	HAH302	HAH303	HAH304	HAH305	HAH306	HAH307	HAH308	HAH309	HAH310	HAH311	HAH312	HAH313	HAH314	HAH315	HAH316	HAH317	HAH318	HAH319	HAH320	HAH321	HAH322	HAH323	HAH324	HAH325	HAH326	HAH327	HAH328	HAH329	HAH330	HAH331	HAH332	HAH333	HAH334	HAH335	HAH336	HAH337	HAH338	HAH339	HAH340	HAH341	HAH342	HAH343	HAH344	HAH345	HAH346	HAH347	HAH348	HAH349	HAH350	HAH351	HAH352	HAH353	HAH354	HAH355	HAH356	HAH357	HAH358	HAH359	HAH360	HAH361	HAH362	HAH363	HAH364	HAH365	HAH366	HAH367	HAH368	HAH369	HAH370	HAH371	HAH372	HAH373	HAH374	HAH375	HAH376	HAH377	HAH378	HAH379	HAH380	HAH381	HAH382	HAH383	HAH384	HAH385	HAH386	HAH387	HAH388	HAH389	HAH390	HAH391	HAH392	HAH393	HAH394	HAH395	HAH396	HAH397	HAH398	HAH399	HAH400	HAH401	HAH402	HAH403	HAH404	HAH405	HAH406	HAH407	HAH408	HAH409	HAH410	HAH411	HAH412	HAH413	HAH414	HAH415	HAH416	HAH417	HAH418	HAH419	HAH420	HAH421	HAH422	HAH423	HAH424	HAH425	HAH426	HAH427	HAH428	HAH429	HAH430	HAH431	HAH432	HAH433	HAH434	HAH435	HAH436	HAH437	HAH438	HAH439	HAH440	HAH441	HAH442	HAH443	HAH444	HAH445	HAH446	HAH447	HAH448	HAH449	HAH450	HAH451	HAH452	HAH453	HAH454	HAH455	HAH456	HAH457	HAH458	HAH459	HAH460	HAH461	HAH462	HAH463	HAH464	HAH465	HAH466	HAH467	HAH468	HAH469	HAH470	HAH471	HAH472	HAH473	HAH474	HAH475	HAH476	HAH477	HAH478	HAH479	HAH480	HAH481	HAH482	HAH483	HAH484	HAH485	HAH486	HAH487	HAH488	HAH489	HAH490	HAH491	HAH492	HAH493	HAH494	HAH495	HAH496	HAH497	HAH498	HAH499	HAH500	HAH501	HAH502	HAH503	HAH504	HAH505	HAH506	HAH507	HAH508	HAH509	HAH510	HAH511	HAH512	HAH513	HAH514	HAH515	HAH516	HAH517	HAH518	HAH519	HAH520	HAH521	HAH522	HAH523	HAH524	HAH525	HAH526	HAH527	HAH528	HAH529	HAH530	HAH531	HAH532	HAH533	HAH534	HAH535	HAH536	HAH537	HAH538	HAH539	HAH540	HAH541	HAH542	HAH543	HAH544	HAH545	HAH546	HAH547	HAH548	HAH549	HAH550	HAH551	HAH552	HAH553	HAH554	HAH555	HAH556	HAH557	HAH558	HAH559	HAH560	HAH561	HAH562	HAH563	HAH564	HAH565	HAH566	HAH567	HAH568	HAH569	HAH570	HAH571	HAH572	HAH573	HAH574	HAH575	HAH576	HAH577	HAH578	HAH579	HAH580	HAH581	HAH582	HAH583	HAH584	HAH585	HAH586	HAH587	HAH588	HAH589	HAH590	HAH591	HAH592	HAH593	HAH594	HAH595	HAH596	HAH597	HAH598	HAH599	HAH600	HAH601	HAH602	HAH603	HAH604	HAH605	HAH606	HAH607	HAH608	HAH609	HAH610	HAH611	HAH612	HAH613	HAH614	HAH615	HAH616	HAH617	HAH618	HAH619	HAH620	HAH621	HAH622	HAH623	HAH624	HAH625	HAH626	HAH627	HAH628	HAH629	HAH630	HAH631	HAH632	HAH633	HAH634	HAH635	HAH636	HAH637	HAH638	HAH639	HAH640	HAH641	HAH642	HAH643	HAH644	HAH645	HAH646	HAH647	HAH648	HAH649	HAH650	HAH651	HAH652	HAH653	HAH654	HAH655	HAH656	HAH657	HAH658	HAH659	HAH660	HAH661	HAH662	HAH663	HAH664	HAH665	HAH666	HAH667	HAH668	HAH669	HAH670	HAH671	HAH672	HAH673	HAH674	HAH675	HAH676	HAH677	HAH678	HAH679	HAH680	HAH681	HAH682	HAH683	HAH684	HAH685	HAH686	HAH687	HAH688	HAH689	HAH690	HAH691	HAH692	HAH693	HAH694	HAH695	HAH696	HAH697	HAH698	HAH699	HAH700	HAH701	HAH702	HAH703	HAH704	HAH705	HAH706	HAH707	HAH708	HAH709	HAH710	HAH711	HAH712	HAH713	HAH714	HAH715	HAH716	HAH717	HAH718	HAH719	HAH720	HAH721	HAH722	HAH723	HAH724	HAH725	HAH726	HAH727	HAH728	HAH729	HAH730	HAH731	HAH732	HAH733	HAH734	HAH735	HAH736	HAH737	HAH738	HAH739	HAH740	HAH741	HAH742	HAH743	HAH744	HAH745	HAH746	HAH747	HAH748	HAH749	HAH750	HAH751	HAH752	HAH753	HAH754	HAH755	HAH756	HAH757	HAH758	HAH759	HAH760	HAH761	HAH762	HAH763	HAH764	HAH765	HAH766	HAH767	HAH768	HAH769	HAH770	HAH771	HAH772	HAH773	HAH774	HAH775	HAH776	HAH777	HAH778	HAH779	HAH780	HAH781	HAH782	HAH783	HAH784	HAH785	HAH786	HAH787	HAH788	HAH789	HAH790	HAH791	HAH792	HAH793	HAH794	HAH795	HAH796	HAH797	HAH798	HAH799	HAH800	HAH801	HAH802	HAH803	HAH804	HAH805	HAH806	HAH807	HAH808	HAH809	HAH810	HAH811	HAH812	HAH813	HAH814	HAH815	HAH816	HAH817	HAH818	HAH819	HAH820	HAH821	HAH822	HAH823	HAH824	HAH825	HAH826	HAH827	HAH828	HAH829	HAH830	HAH831	HAH832	HAH833	HAH834	HAH835	HAH836	HAH837	HAH838	HAH839	HAH840	HAH841	HAH842	HAH843	HAH844	HAH845	HAH846	HAH847	HAH848	HAH849	HAH850	HAH851	HAH852	HAH853	HAH854	HAH855	HAH856	HAH857	HAH858	HAH859	HAH860	HAH861	HAH862	HAH863	HAH864	HAH865	HAH866	HAH867	HAH868	HAH869	HAH870	HAH871	HAH872	HAH873	HAH874	HAH875	HAH876	HAH877	HAH878	HAH879	HAH880	HAH881	HAH882	HAH883	HAH884	HAH885	HAH886	HAH887	HAH888	HAH889	HAH890	HAH891	HAH892	HAH893	HAH894	HAH895	HAH896	HAH897	HAH898	HAH899	HAH900	HAH901	HAH902	HAH903	HAH904	HAH905	HAH906	HAH907	HAH908	HAH909	HAH910	HAH911	HAH912	HAH913	HAH914	HAH915	HAH916	HAH917	HAH918	HAH919	HAH920	HAH921	HAH922	HAH923	HAH924	HAH925	HAH926	HAH927	HAH928	HAH929	HAH930	HAH931	HAH932	HAH933	HAH934	HAH935	HAH936	HAH937	HAH938	HAH939	HAH940	HAH941	HAH942	HAH943	HAH944	HAH945	HAH946	HAH947	HAH948	HAH949	HAH950	HAH951	HAH952	HAH953	HAH954	HAH955	HAH956	HAH957	HAH958	HAH959	HAH960	HAH961	HAH962	HAH963	HAH964	HAH965	HAH966	HAH967	HAH968	HAH969	HAH970	HAH971	HAH972	HAH973	HAH974	HAH975	HAH976	HAH977	HAH978	HAH979	HAH980	HAH981	HAH982	HAH983	HAH984	HAH985	HAH986	HAH987	HAH988	HAH989	HAH990	HAH991	HAH992	HAH993	HAH994	HAH995	HAH996	HAH997	HAH998	HAH999	HAH1000

MC46423: Chain of Custody

Page 1 of 3

SGS

25 of 811
ACCUTEST
MC46423

EXECUTIVE NARRATIVE

SDG No: **MC46423** Laboratory: **Accutest, Massachusetts**
Analysis: **MADEP VPH** Number of Samples: **7**
Location: **BMSMC, Building 5 Area**
Humacao, PR

SUMMARY: Seven (7) samples were analyzed for Volatiles TPHC Ranges by method MADEP VPH. Samples were validated following the METHOD FOR THE DETERMINATION OF VOLATILE PETROLEUM HYDROCARBONS (VPH) quality control criteria, Massachusetts Department of Environmental Protection, Revision 1.1 (2004). Also the general validation guidelines promulgated by the USEPA Hazardous Wastes Support Section. The QC criteria and data validation actions listed on the data review worksheets are from the primary guidance document, unless otherwise noted.

Results are valid and can be used for decision making purposes.

Critical issues: **None**
Major: **None**
Minor: **None**

Critical findings: **None**
Major findings: **None**
Minor findings: **None**

COMMENTS: Results are valid and can be used for decision making purposes.

Reviewers Name: **Rafael Infante**
Chemist License 1888

Signature:

A handwritten signature in blue ink, reading "Rafael Infante", is written over a horizontal line.

Date: **July 16, 2016**

SAMPLE ORGANIC DATA SAMPLE SUMMARY

Sample ID: MC46423-1
Sample location: BMSMC Building 5 Area
Sampling date: 6/13/2016
Matrix: AQ - Equipment Blank

METHOD: MADEP VPH

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Ç5 - C8 Aliphatics (Unadj.)	30.9	ug/L	1	J	UJ	Yes
Ç9 - C12 Aliphatics (Unadj.)	72.3	ug/L	1	-	-	Yes
Ç9 - C10 Aromatics (Unadj.)	50.0	ug/L	1	-	-	Yes
Ç5 - C8 Aliphatics	50	ug/L	1	-	U	Yes
Ç9 - C12 Aliphatics	50	ug/L	1	-	U	Yes

Sample ID: MC46423-2
Sample location: BMSMC Building 5 Area
Sampling date: 6/13/2016
Matrix: Groundwater

METHOD: MADEP VPH

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Ç5 - C8 Aliphatics (Unadj.)	50	ug/L	1	-	U	Yes
Ç9 - C12 Aliphatics (Unadj.)	50	ug/L	1	-	U	Yes
Ç9 - C10 Aromatics (Unadj.)	50	ug/L	1	-	U	Yes
Ç5 - C8 Aliphatics	50	ug/L	1	-	U	Yes
Ç9 - C12 Aliphatics	50	ug/L	1	-	U	Yes

Sample ID: MC46423-3
Sample location: BMSMC Building 5 Area
Sampling date: 6/15/2016
Matrix: Groundwater

METHOD: MADEP VPH

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Ç5 - C8 Aliphatics (Unadj.)	167	ug/L	1	-	-	Yes
Ç9 - C12 Aliphatics (Unadj.)	63100	ug/L	100	-	-	Yes
Ç9 - C10 Aromatics (Unadj.)	112	ug/L	1	-	-	Yes
Ç5 - C8 Aliphatics	67.1	ug/L	1	-	-	Yes
Ç9 - C12 Aliphatics	753	ug/L	1	-	-	Yes

Sample ID: MC46423-4
Sample location: BMSMC Building 5 Area
Sampling date: 6/15/2016
Matrix: Groundwater

METHOD: MADEP VPH

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Ç5 - C8 Aliphatics (Unadj.)	36.4	ug/L	1	J	UJ	Yes
Ç9 - C12 Aliphatics (Unadj.)	55.4	ug/L	1	-	-	Yes
Ç9 - C10 Aromatics (Unadj.)	50	ug/L	1	-	U	Yes
Ç5 - C8 Aliphatics	27.7	ug/L	1	J	UJ	Yes
Ç9 - C12 Aliphatics	50	ug/L	1	-	U	Yes

Sample ID: MC46423-5
Sample location: BMSMC Building 5 Area
Sampling date: 6/15/2016
Matrix: Groundwater

METHOD: MADEP VPH

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Ç5 - C8 Aliphatics (Unadj.)	50	ug/L	1	-	U	Yes
Ç9 - C12 Aliphatics (Unadj.)	50	ug/L	1	-	U	Yes
Ç9 - C10 Aromatics (Unadj.)	50	ug/L	1	-	U	Yes
Ç5 - C8 Aliphatics	50	ug/L	1	-	U	Yes
Ç9 - C12 Aliphatics	50	ug/L	1	-	U	Yes

Sample ID: MC46423-6
Sample location: BMSMC Building 5 Area
Sampling date: 6/16/2016
Matrix: Groundwater

METHOD: MADEP VPH

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Ç5 - C8 Aliphatics (Unadj.)	27.2	ug/L	1	J	UJ	Yes
Ç9 - C12 Aliphatics (Unadj.)	19500	ug/L	1	-	-	Yes
Ç9 - C10 Aromatics (Unadj.)	100	ug/L	1	-	-	Yes
Ç5 - C8 Aliphatics	26.3	ug/L	1	J	UJ	Yes
Ç9 - C12 Aliphatics	65.3	ug/L	1	-	-	Yes

Sample ID: MC46423-7
Sample location: BMSMC Building 5 Area
Sampling date: 6/16/2016
Matrix: Groundwater

METHOD: MADEP VPH

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Ç5 - C8 Aliphatics (Unadj.)	50	ug/L	1	-	UJ	Yes
Ç9 - C12 Aliphatics (Unadj.)	50	ug/L	1	-	U	Yes
Ç9 - C10 Aromatics (Unadj.)	50	ug/L	1	-	U	Yes
Ç5 - C8 Aliphatics	50	ug/L	1	-	U	Yes
Ç9 - C12 Aliphatics	50	ug/L	1	-	U	Yes

DATA REVIEW WORKSHEETS

Type of validation Full: ☒ _____
 Limited: _____
 Project Number: MC46423
 Date: 06/13-16/2016
 Shipping date: 06/16/2016
 EPA Region: 2

REVIEW OF VOLATILE PETROLEUM HYDROCARBON (VPHs) PACKAGE

The following guidelines for evaluating volatile organics were created to delineate required validation actions. This document will assist the reviewer in using professional judgment to make more informed decision and in better serving the needs of the data users. The sample results were assessed according to the data validation guidance documents in the following order of precedence METHOD FOR THE DETERMINATION OF VOLATILE PETROLEUM HYDROCARBONS (VPH), Massachusetts Department of Environmental Protection, Revision 1.1 (2004). Also the general validation guidelines promulgated by the USEPA Hazardous Wastes Support Section. The QC criteria and data validation actions listed on the data review worksheets are from the primary guidance document, unless otherwise noted.

The hardcopied (laboratory name) Accutest Laboratories data package received has been reviewed and the quality control and performance data summarized. The data review for SVOCs included:

Lab. Project/SDG No.: MC46423 Sample matrix: Groundwater
 No. of Samples: 7
 Field blank No.: _____
 Equipment blank No.: _____
 Trip blank No.: _____
 Field duplicate No.: _____

<input checked="" type="checkbox"/> Data Completeness	<input checked="" type="checkbox"/> Laboratory Control Spikes
<input checked="" type="checkbox"/> Holding Times	<input checked="" type="checkbox"/> Field Duplicates
<input type="checkbox"/> N/A GC/MS Tuning	<input checked="" type="checkbox"/> Calibrations
<input type="checkbox"/> N/A Internal Standard Performance	<input checked="" type="checkbox"/> Compound Identifications
<input checked="" type="checkbox"/> Blanks	<input checked="" type="checkbox"/> Compound Quantitation
<input checked="" type="checkbox"/> Surrogate Recoveries	<input checked="" type="checkbox"/> Quantitation Limits
<input checked="" type="checkbox"/> Matrix Spike/Matrix Spike Duplicate	

Overall Comments: Volatiles by GC by Method MADEP_VPH_REV_1.1
(C5_to_C12_Aliphatics; C9_to_C10_Aromatics)

Definition of Qualifiers:

J- Estimated results
 U- Compound not detected
 R- Rejected data
 UJ- Estimated nondetect

Reviewer: Rafael Infante
 Date: 07/16/2016

DATA REVIEW WORKSHEETS

All criteria were met x
Criteria were not met and/or see below

I. DATA COMPLETNESS

A. Data Package:

MISSING INFORMATION

DATE LAB. CONTACTED

DATE RECEIVED

[illegible]

B. Other

Discrepancies:

[The page contains faint horizontal lines, suggesting it was part of a lined notebook or document.]

DATA REVIEW WORKSHEETS

All criteria were met X
Criteria were not met and/or see below

HOLDING TIMES

The objective of this parameter is to ascertain the validity of the results based on the holding time of the sample from time of collection to the time of extraction, and subsequently from the time of extraction to the time of analysis.

Complete table for all samples and note the analysis and/or preservation not within criteria

SAMPLE ID	DATE SAMPLED	DATE EXTRACTED	DATE ANALYZED	ACTION
Samples analyzed within method recommended holding time				

Criteria

Preservation:

Samples analyzed with ambient purge temperature: Samples must be acidified to a pH of 2.0 or less at the time of collection.

Samples analyzed with heated purge temperature: Samples must be treated to a pH of 11.0 or greater at the time of collection.

Methanol preservation of soil/sediment samples is mandatory. Methanol (purge-and-trap grade) must be added to the sample vial before or immediately after sample collection. In lieu of the in-field preservation of samples with methanol, soil samples may be obtained in specially-designed air tight sampling devices, provided that the samples are extruded and preserved in methanol within 48 hours of collection.

Holding times:

Aqueous samples using ambient or heated purge - analyze within 14 days.

Soil/sediment samples - analysis within 28 days.

Cooler temperature (Criteria: 4 ± 2 °C): 2.5°C

Actions: Qualify positive results/non-detects as follows:

If holding times are exceeded, estimate positive results (J) and nondetects (UJ).

If holding times are grossly exceeded, use professional judgment to qualify data. The data reviewer may choose to estimate positive results (J) and rejects nondetects (R).

If samples were not at the proper temperature ($> 10^{\circ}\text{C}$) or improperly preserved, use professional judgment to qualify the results.

DATA REVIEW WORKSHEETS

All criteria were met X
Criteria were not met and/or see below

CALIBRATIONS VERIFICATION

Compliance requirements for satisfactory instrument calibration are established to ensure that the instrument is capable of producing and maintaining acceptable quantitative data.

Date of initial calibration: 06/02/16

Dates of initial calibration verification: 06/02/16

Instrument ID numbers: GCWX

Matrix/Level: AQUEOUS/MEDIUM

DATE	LAB FILE ID#	ANALYTE	CRITERIA OUT RFs, %RSD, %D, r	SAMPLES AFFECTED
Initial and initial calibration verification meet method specific requirements				

Criteria- ICAL

- Five point calibration curve.
- The percent relative standard deviation (%RSD) of the calibration factor must be equal to or less than 25% over the working range for the analyte of interest. When this condition is met, linearity through the origin may be assumed, and the average calibration factor is used in lieu of a calibration curve.
- A collective calibration factor must also be established for each hydrocarbon range of interest. Calculate the collective CFs for C5-C8 Aliphatic Hydrocarbons and C9-C12 Aliphatic Hydrocarbons using the FID chromatogram. Calculate the collective CF for the C9-C10 Aromatic Hydrocarbons using the PID chromatogram. Tabulate the summation of the peak areas of all components in that fraction against the total concentration injected. The %RSD of the calibration factor must be equal to or less than 25% over the working range for the hydrocarbon range of interest.

Criteria- CCAL

- At a minimum, the working calibration factor must be verified on each working day, after every 20 samples, and at the end of the analytical sequence by the injection of a mid-level continuing calibration standard to verify instrument performance and linearity.
- If the percent difference (%D) for any analyte varies from the predicted response by more than $\pm 25\%$, a new five-point calibration must be performed for that analyte. Greater percent differences are permissible for n-nonane. If the %D for n-nonane is greater than 30, note the nonconformance in the case narrative. It

DATA REVIEW WORKSHEETS

should be noted that the %Ds are calculated when CFs are used for the initial calibration and percent drifts are calculated when calibration curves using linear regression are used for the initial calibration.

Actions:

If %RSD > 25% for target compounds or a correlation coefficient < 0.99, estimate positive results (J) and use professional judgment to qualify nondetects.

If % D > 25% (> 30 for nonane), estimate positive results (J) and nondetects (UJ).

CALIBRATIONS VERIFICATION

Compliance requirements for satisfactory instrument calibration are established to ensure that the instrument is capable of producing and maintaining acceptable quantitative data.

Date of initial calibration: _____ 06/09/16 _____

Dates of continuing calibration verification: _____ 06/20/16 _____

Dates of final calibration verification: _____ 06/21/16 _____

Instrument ID numbers: _____ GCWX _____

Matrix/Level: _____ AQUEOUS/MEDIUM _____

DATE	LAB FILE ID#	ANALYTE	CRITERIA OUT RFs, %RSD, %D, r	SAMPLES AFFECTED
Continuing and final calibration verification meet method specific requirements				

A separate worksheet should be filled for each initial curve

DATA REVIEW WORKSHEETS

All criteria were met X
Criteria were not met and/or see below

V A. BLANK ANALYSIS RESULTS (Sections 1 & 2)

The assessment of the blank analysis results is to determine the existence and magnitude of contamination problems. The criteria for evaluation of blanks apply only to blanks associated with the samples, including trip, equipment, and laboratory blanks. If problems with any blanks exist, all data associated with the case must be carefully evaluated to determine whether or not there is an inherent variability in the data for the case, or if the problem is an isolated occurrence not affecting other data. A Laboratory Method Blank must be run after samples suspected of being highly contaminated to determine if sample carryover has occurred.

List the contamination in the blanks below. High and low levels blanks must be treated separately.

Laboratory blanks

DATE ANALYZED	LAB ID	LEVEL/MATRIX	COMPOUND	CONCENTRATION UNITS
---------------	--------	--------------	----------	---------------------

 METHOD BLANKS MEET THE METHOD SPECIFIC CRITERIA

Field/Trip/Equipment

A methanol trip blank or acidified reagent water trip blank should continually accompany each soil/sediment sample or water sample batch, respectively, during sampling, storage, and analysis.

DATE ANALYZED	LAB ID	LEVEL/MATRIX	COMPOUND	CONCENTRATION UNITS
---------------	--------	--------------	----------	---------------------

 NO TRIP/FIELD/EQUIPMENT BLANKS ASSOCIATED WITH THIS DATA PACKAGE

DATA REVIEW WORKSHEETS

All criteria were met X
Criteria were not met and/or see below

V B. BLANK ANALYSIS RESULTS (Section 3)

Blank Actions

The ALs for samples which have been diluted should be corrected for the sample dilution factor and/or % moisture, where applicable. Peaks must not be detected above the Reporting Limit within the retention time window of any analyte of interest. The hydrocarbon ranges must not be detected at a concentration greater than 10% of the most stringent MCP cleanup standard. Specific actions area as follows:

If the concentration is $<$ sample quantitation limit (SQL) and $<$ AL, report the compound as not detected (U) at the SQL.

If the concentration is \geq SQL but $<$ AL, report the compound as not detected (U) at the reported concentration.

If the concentration is $>$ AL, report the concentration unqualified.

DATA REVIEW WORKSHEETS

All criteria were met X
Criteria were not met and/or see below

SURROGATE SPIKE RECOVERIES

Laboratory performance of individual samples is established by evaluation of surrogate spike recoveries. All samples are spiked with surrogate compounds prior to sample analysis. The accuracy of the analysis is measured by the surrogate percent recovery. Since the effects of the sample matrix are frequently outside the control of the laboratory and may present relatively unique problems, the validation of data is frequently subjective and demands analytical experience and professional judgment.

List the percent recoveries (%Rs) which do not meet the criteria for surrogate recovery.

Matrix: solid/aqueous

SAMPLE ID	SURROGATE COMPOUND	ACTION
	2,3,4-Trifluorotoluene	
<u> SURROGATE_STANDARD_RECOVERIES_WITHIN_LABORATORY_CONTROL </u>		
<u> LIMITS </u>		

QC Limits* (Aqueous)

 LL to UL 70 to 130 to to

QC Limits* (Solid)

 LL to UL to to to

It is recommended that surrogate standard recoveries be monitored and documented on a continuing basis. At a minimum, when surrogate recovery from a sample, blank, or QC sample is less than 70% or more than 130%, check calculations to locate possible errors, check the fortifying standard solution for degradation, and check changes in instrument performance.

If the cause cannot be determined, reanalyze the sample unless one of the following exceptions applies:

- (1) Obvious interference is present on the chromatogram (e.g., unresolved complex mixture);
- (2) Percent moisture of associated soil/sediment sample is >25% and surrogate recovery is >10%; or
- (3) The surrogate exhibits high recovery and associated target analytes or hydrocarbon ranges are not detected in sample.

If a sample with a surrogate recovery outside of the acceptable range is not reanalyzed based on any of these aforementioned exceptions, this information must be noted on the data report form and discussed in the Executive Report. Analysis of the sample on dilution may diminish matrix-related surrogate recovery problems. This approach can be used as long as the reporting limits to evaluate applicable MCP standards can still be achieved with the dilution. If not, reanalysis without dilution must be performed.

DATA REVIEW WORKSHEETS

All criteria were met X
 Criteria were not met and/or see below

VII. A MATRIX SPIKE/MATRIX SPIKE DUPLICATE (MS/MSD)

This data is generated to determine long term precision and accuracy in the analytical method for various matrices. This data alone cannot be used to evaluate the precision and accuracy of individual samples.

At the request of the data user, and in consideration of sample matrices and data quality objectives, matrix spikes and matrix duplicates may be analyzed with every batch of 20 samples or less per matrix.

- **Matrix duplicate** - Matrix duplicates are prepared by analyzing one sample in duplicate. The purpose of the matrix duplicates is to determine the homogeneity of the sample matrix as well as analytical precision. The RPD of detected results in the matrix duplicate samples must not exceed 50 when the results are greater than 5x the reporting limit.
- The desired spiking level is 50% of the highest calibration standard. However, the total concentration in the MS (including the MS and native concentration in the unspiked sample) should not exceed 75% of the highest calibration standard in order for a proper evaluation to be performed. The purpose of the matrix spike is to determine whether the sample matrix contributes bias to the analytical results. The corrected concentrations of each analyte within the matrix spiking solution must be within 70 - 130% of the true value. Lower recoveries of n-nonane are permissible (if included in the calibration of the C9-C12 aliphatic range), but must be noted in the narrative if <30%.

MS/MSD Recoveries and Precision Criteria

Sample ID: MC46423-7_MS/MSD Matrix/Level: Groundwater

List the %Rs, RPD of the compounds which do not meet the QC criteria.

MS OR MSD	COMPOUND	% R	RPD	QC LIMITS	ACTION

Note: MS/MSD % recoveries and RPD within laboratory control limits.

DATA REVIEW WORKSHEETS

All criteria were met X
 Criteria were not met and/or see below

No action is taken on MS/MSD results alone to qualify the entire case. However, used informed professional judgment, the data reviewer may use the MS/MSD results in conjunction with other QC criteria and determine the need for some qualification of the data. In those instances where it can be determined that the results of the MS/MSD affect only the sample spiked, the qualification should be limited to this sample alone. However, it may be determined through the MS/MSD results that the laboratory is having a systematic problem in the analysis of one or more analytes, which affects the associated samples.

2. MS/MSD – Unspiked Compounds

List the concentrations of the unspiked compounds and determine the % RSDs of these compounds in the unspiked sample, matrix spike, and matrix spike duplicate.

COMPOUND	CONCENTRATION			%RPD	ACTION
	SAMPLE	MS	MSD		

Criteria: None specified, use %RSD \leq 50 as professional judgment.

Actions:

If the % RSD > 50, qualify the results in the spiked sample as estimate (J).

If the % RSD is not calculable (NC) due to nondetect value in the sample, MS, and/or MSD, use professional judgment to qualify sample data.

A separate worksheet should be used for each MS/MSD pair.

DATA REVIEW WORKSHEETS

All criteria were met __X__
Criteria were not met and/or see below _____

VIII. LABORATORY CONTROL SAMPLE (LCS/LCSD) ANALYSIS

This data is generated to determine accuracy of the analytical method for various matrices.

1. LCS Recoveries Criteria

List the %R of compounds which do not meet the criteria

LCS ID	COMPOUND	% R	QC LIMIT	ACTION
--------	----------	-----	----------	--------

__LCS_RECOVERY_WITHIN_LABORATORY_CONTROL_LIMITS__

Criteria:

- * Refer to QAPP for specific criteria.
- * The spike recovery must be between 70% and 130%. Lower recoveries of n-nonane are permissible (if included in the calibration of the C9-C12 aliphatic range). If the recovery of n-nonane is <30%, note the nonconformance in the executive narrative.

Actions:

Actions on LCS recovery should be based on both the number of compounds that are outside the %R criteria and the magnitude of the exceedance of the criteria.

If the %R of the analyte is > UL, qualify all positive results (j) for the affected analyte in the associated samples and accept nondetects.

If the %R of the analyte is < LL, qualify all positive results (j) and reject (R) nondetects for the affected analyte in the associated samples.

If more than half the compounds in the LCS are not within the required recovery criteria, qualify all positive results as (J) and reject nondetects (R) for all target analyte(s) in the associated samples.

2. Frequency Criteria:

Where LCS analyzed at the required frequency and for each matrix (1 per 20 samples per matrix)? Yes or No.

If no, the data may be affected. Use professional judgment to determine the severity of the effect and qualify data accordingly. Discuss any actions below and list the samples affected. Discuss the actions below:

DATA REVIEW WORKSHEETS

All criteria were met N/A
Criteria were not met and/or see below

IX. FIELD/LABORATORY DUPLICATE PRECISION

Sample IDs: -

Matrix: -

Field/laboratory duplicate samples may be taken and analyzed as an indication of overall precision. These analyses measure both field and lab precision; therefore, the results may have more variability than laboratory duplicates which measures only laboratory performance. It is also expected that soil duplicate results will have a greater variance than water matrices due to difficulties associated with collecting identical field duplicate samples.

COMPOUND	SQL	SAMPLE CONC.	DUPLICATE CONC.	RPD	ACTION
No field/laboratory duplicate analyzed with this data package. MS/MSD % recovery RPD used to assess accuracy. RPD within laboratory and validation guidance document criteria (+ 50 %) for analytes detected above reporting limits.					

Criteria:

The project QAPP should be reviewed for project-specific information.
RPD \pm 30% for aqueous samples, RPD \pm 50 % for solid samples if results are \geq SQL.
If both samples and duplicate are < 5 SQL, the RPD criteria is doubled.

SQL = soil quantitation limit

Actions:

If both the sample and the duplicate results are nondetects (ND), the RPD is not calculable (NC). No action is needed.

Qualify as estimated positive results (J) and nondetects (UJ) for the compound that exceeded the above criteria.

If one sample result is not detected and the other is ≥ 5 x the SQL qualify (J/UJ).

Note: If SQLs for the sample and duplicate are significantly different, use professional judgment to determine if qualification is appropriate.

If one sample value is not detected and the other is < 5 x the SQL, use professional judgment to determine if qualification is appropriate.

DATA REVIEW WORKSHEETS

All criteria were met X
Criteria were not met and/or see below

XI. COMPOUND IDENTIFICATION

The compound identification evaluation is to verify that the laboratory correctly identified target analytes as well as tentatively identified compounds (TICs).

1. Verify that the target analytes were within the retention time windows.
 - Retention time windows must be re-established for each Target VPH Analyte each time a new GC column is installed, and must be verified and/or adjusted on a daily basis.
 - Coelution of the m- and p- xylene isomers is permissible.
 - All surrogates must be adequately resolved from individual Target Analytes included in the VPH Component Standard.
 - For the purposes of this method, adequate resolution is assumed to be achieved if the height of the valley between two peaks is less than 25% of the average height of the two peaks.
 - The n-pentane (C5) and MtBE peaks must be adequately resolved from any solvent front that may be present on the FID and PID chromatograms, respectively.

Note: Target analytes were within the retention time window.

2. If target analytes and/or TICs were not correctly identified, request that the laboratory resubmit the corrected data.

DATA REVIEW WORKSHEETS

All criteria were met X
Criteria were not met and/or see below

XII. QUANTITATION LIMITS AND SAMPLE RESULTS

The sample quantitation evaluation is to verify laboratory quantitation results.

1. In the space below, please show a minimum of one sample calculation:

MC46423-3 VPH (C5 – C7 Aliphatics) RF = 2.366×10^4

FID

$$[] = (24176) / (2.366 \times 10^4)$$

$$[] = 1.02 \text{ ppb} \quad \text{Ok}$$

MC46423-1 VPH (C9 – C10 Aromatics) RF = 1.264×10^4

PID

$$[] = (1412289) / (1.264 \times 10^4)$$

$$[] = 111.7 \text{ ppb} \quad \text{Ok}$$

2. If requested, verify that the results were above the laboratory method detection limit (MDLs).
3. If dilutions performed, were the SQLs elevated accordingly by the laboratory? List the affected samples and dilution factor in the table below.

SAMPLE ID	DILUTION FACTOR	REASON FOR DILUTION
MC46423-3	100 X	C9 – C12 aliphatic hydrocarbon range over calibration range
MC46423-6	50 X	

If dilution was not performed and the results were above the concentration range, estimate results (J) for the affected compounds. List the affected samples/compounds:

EXECUTIVE NARRATIVE

SDG No: **MC46423** Laboratory: **Accutest, Massachusetts**
Analysis: **MADEP EPH** Number of Samples: **7**
Location: **BMSMC, Building 5 Area**
Humacao, PR

SUMMARY: Seven (7) samples were analyzed for Extractable Petroleum Hydrocarbons TPHC Ranges by method MADEP EPH. Samples were validated following the METHOD FOR THE DETERMINATION OF EXTRACTABLE PETROLEUM HYDROCARBONS (EPH) quality control criteria, Massachusetts Department of Environmental Protection, Revision 1.1 (2004). Also the general validation guidelines promulgated by the USEPA Hazardous Wastes Support Section. The QC criteria and data validation actions listed on the data review worksheets are from the primary guidance document, unless otherwise noted.

Results are valid and can be used for decision making purposes.

Critical issues: **None**
Major: **None**
Minor: **None**

Critical findings: **None**
Major findings: **None**
Minor findings: **None**

COMMENTS: Results are valid and can be used for decision making purposes.

Reviewers Name: **Rafael Infante**
Chemist License 1888

Signature:

A handwritten signature in blue ink, reading "Rafael Infante", is written over a horizontal line.

Date: **July 16, 2016**

SAMPLE ORGANIC DATA SAMPLE SUMMARY

Sample ID: MC46423-1A
Sample location: BMSMC Building 5 Area
Sampling date: 6/13/2016
Matrix: Groundwater

METHOD: MADEP EPH

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Ç11 - C22 Aromatics (Unadj.)	33.2	ug/L	1	JB	UJ	Yes
Ç9 - C18 Aliphatics	20.6	ug/L	1	J	UJ	Yes
Ç19 - C36 Aliphatics	32.9	ug/L	1	J	UJ	Yes
Ç11 - C22 Aromatics	32.6	ug/L	1	JB	UJ	Yes

Sample ID: MC46423-2A
Sample location: BMSMC Building 5 Area
Sampling date: 6/13/2016
Matrix: Groundwater

METHOD: MADEP EPH

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Ç11 - C22 Aromatics (Unadj.)	34.2	ug/L	1	JB	UJ	Yes
Ç9 - C18 Aliphatics	22.7	ug/L	1	J	UJ	Yes
Ç19 - C36 Aliphatics	40.9	ug/L	1	J	UJ	Yes
Ç11 - C22 Aromatics	34.2	ug/L	1	JB	UJ	Yes

Sample ID: MC46423-3A
Sample location: BMSMC Building 5 Area
Sampling date: 6/15/2016
Matrix: Groundwater

METHOD: MADEP EPH

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Ç11 - C22 Aromatics (Unadj.)	36.4	ug/L	1	JB	UJ	Yes
Ç9 - C18 Aliphatics	25.8	ug/L	1	J	UJ	Yes
Ç19 - C36 Aliphatics	96.7	ug/L	1	J	UJ	Yes
Ç11 - C22 Aromatics	36.4	ug/L	1	JB	UJ	Yes

Sample ID: MC46423-4A
Sample location: BMSMC Building 5 Area
Sampling date: 6/15/2016
Matrix: Groundwater

METHOD: MADEP EPH

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Ç11 - C22 Aromatics (Unadj.)	60.6	ug/L	1	JB	UJ	Yes
Ç9 - C18 Aliphatics	22.7	ug/L	1	J	UJ	Yes
Ç19 - C36 Aliphatics	46.6	ug/L	1	J	UJ	Yes
Ç11 - C22 Aromatics	32.4	ug/L	1	JB	UJ	Yes

Sample ID: MC46423-5A
Sample location: BMSMC Building 5 Area
Sampling date: 6/15/2016
Matrix: Groundwater

METHOD: MADEP EPH

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Ç11 - C22 Aromatics (Unadj.)	110	ug/L	1	-	U	Yes
Ç9 - C18 Aliphatics	20.8	ug/L	1	J	UJ	Yes
Ç19 - C36 Aliphatics	33.7	ug/L	1	J	UJ	Yes
Ç11 - C22 Aromatics	110	ug/L	1	-	U	Yes

Sample ID: MC46423-6A
Sample location: BMSMC Building 5 Area
Sampling date: 6/16/2016
Matrix: Groundwater

METHOD: MADEP EPH

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Ç11 - C22 Aromatics (Unadj.)	85.3	ug/L	1	JB	UJ	Yes
Ç9 - C18 Aliphatics	72.7	ug/L	1	J	UJ	Yes
Ç19 - C36 Aliphatics	37.3	ug/L	1	J	UJ	Yes
Ç11 - C22 Aromatics	78.4	ug/L	1	JB	UJ	Yes

Sample ID: MC46423-7A
Sample location: BMSMC Building 5 Area
Sampling date: 6/16/2016
Matrix: Groundwater

METHOD: MADEP EPH

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Ç11 - C22 Aromatics (Unadj.)	39.1	ug/L	1	JB	UJ	Yes
Ç9 - C18 Aliphatics	19.9	ug/L	1	J	UJ	Yes
Ç19 - C36 Aliphatics	59.6	ug/L	1	J	UJ	Yes
Ç11 - C22 Aromatics	39.1	ug/L	1	JB	UJ	Yes

DATA REVIEW WORKSHEETS

Type of validation Full: ☒ Limited: ☐ Project Number: MC46423
Date: 06/13-16/2016
Shipping date: 06/16/2016
EPA Region: 2

REVIEW OF EXTRACTABLE PETROLEUM HYDROCARBON (EPHs) PACKAGE

The following guidelines for evaluating volatile organics were created to delineate required validation actions. This document will assist the reviewer in using professional judgment to make more informed decision and in better serving the needs of the data users. The sample results were assessed according to the data validation guidance documents in the following order of precedence METHOD FOR THE DETERMINATION OF EXTRACTABLE PETROLEUM HYDROCARBONS (VPH), Massachusetts Department of Environmental Protection, Revision 1.1 (2004). Also the general validation guidelines promulgated by the USEPA Hazardous Wastes Support Section. The QC criteria and data validation actions listed on the data review worksheets are from the primary guidance document, unless otherwise noted.

The hardcopied (laboratory name) Accutest Laboratories data package received has been reviewed and the quality control and performance data summarized. The data review for SVOCs included:

Lab. Project/SDG No.: MC46423 Sample matrix: Groundwater/Soil
No. of Samples: 7
Field blank No.: -
Equipment blank No.: -
Trip blank No.: -
Field duplicate No.: -

<input checked="" type="checkbox"/> Data Completeness	<input checked="" type="checkbox"/> Laboratory Control Spikes
<input checked="" type="checkbox"/> Holding Times	<input checked="" type="checkbox"/> Field Duplicates
<input type="checkbox"/> GC/MS Tuning	<input checked="" type="checkbox"/> Calibrations
<input type="checkbox"/> Internal Standard Performance	<input checked="" type="checkbox"/> Compound Identifications
<input checked="" type="checkbox"/> Blanks	<input checked="" type="checkbox"/> Compound Quantitation
<input checked="" type="checkbox"/> Surrogate Recoveries	<input checked="" type="checkbox"/> Quantitation Limits
<input checked="" type="checkbox"/> Matrix Spike/Matrix Spike Duplicate	

Overall Extractable Petroleum Hydrocarbons by GC by Method MADEP EPH_REV_1.1 Comments:
(C9_to_C36_Aliphatics; C11_to_C22_Aromatics)

Definition of Qualifiers:

J- Estimated results
U- Compound not detected
R- Rejected data
UJ- Estimated nondetected

Reviewer: Rafael Difant
Date: 07/16/2016

DATA REVIEW WORKSHEETS

All criteria were met x
Criteria were not met and/or see below

I. DATA COMPLETNESS

A. Data Package:

MISSING INFORMATION

DATE LAB. CONTACTEDDATE RECEIVED[illegible]

B. Other

Discrepancies:

[The page contains faint horizontal lines, suggesting it was part of a lined notebook or document.]

DATA REVIEW WORKSHEETS

All criteria were met X
Criteria were not met and/or see below

HOLDING TIMES

The objective of this parameter is to ascertain the validity of the results based on the holding time of the sample from time of collection to the time of extraction, and subsequently from the time of extraction to the time of analysis.

Complete table for all samples and note the analysis and/or preservation not within criteria

SAMPLE ID	DATE SAMPLED	DATE EXTRACTED	DATE ANALYZED	ACTION
Samples extracted and analyzed within method recommended holding time				

Criteria

Preservation:

Aqueous samples must be acidified to a pH of 2.0 or less at the time of collection.

Soil samples must be cooled at 4 ± 2 °C immediately after collection.

Holding times:

Samples must be extracted within 14 days of collection, and analyzed within 40 days of extraction.

Cooler temperature (Criteria: 4 ± 2 °C): 2.5°C

Actions: Qualify positive results/nondetects as follows:

If holding times are exceeded, estimate positive results (J) and nondetects (UJ).

If holding times are grossly exceeded, use professional judgment to qualify data. The data reviewer may choose to estimate positive results (J) and rejects nondetects (R).

If samples were not at the proper temperature ($> 10^{\circ}\text{C}$) or improperly preserved, use professional judgment to qualify the results.

DATA REVIEW WORKSHEETS

All criteria were met X
 Criteria were not met and/or see below

CALIBRATIONS VERIFICATION

Compliance requirements for satisfactory instrument calibration are established to ensure that the instrument is capable of producing and maintaining acceptable quantitative data.

Date of initial calibration: 06/22/16

Dates of initial calibration verification: 06/22/13

Instrument ID numbers: GCDE

Matrix/Level: AQUEOUS/MEDIUM

DATE	LAB FILE ID#	ANALYTE	CRITERIA OUT RFs, %RSD, %D, r	SAMPLES AFFECTED
Initial and continuing calibration meet method specific requirements				

Criteria- ICAL

- Five point calibration curve.
- The percent relative standard deviation (%RSD) of the calibration factor must be equal to or less than 25% over the working range for the analyte of interest. When this condition is met, linearity through the origin may be assumed, and the average calibration factor is used in lieu of a calibration curve.
- A collective calibration factor must also be established for each hydrocarbon range of interest. Calculate the collective CFs for C9-C18 Aliphatic Hydrocarbons, C19-C36 Aliphatic Hydrocarbons, and C11-C22 Aromatic Hydrocarbons using the FID chromatogram. Tabulate the summation of the peak areas of all components in that fraction against the total concentration injected. The %RSD of the calibration factor must be equal to or less than 25% over the working range for the hydrocarbon range of interest.
 - The area for the surrogates must be subtracted from the area summation of the range in which they elute.
 - The areas associated with naphthalene and 2-methylnaphthalene in the aliphatic range standard must be subtracted from the uncorrected collective C9-C18 Aliphatic Hydrocarbon range area prior to calculating the CF.

Criteria- CCAL

- At a minimum, the working calibration factor must be verified on each working day, after every 20 samples or every 24 hours (whichever is more frequent), and

DATA REVIEW WORKSHEETS

at the end of the analytical sequence by the injection of a mid-level continuing calibration standard to verify instrument performance and linearity.

- If the percent difference (%D) for any analyte varies from the predicted response by more than $\pm 25\%$, a new five-point calibration must be performed for that analyte. Greater percent differences are permissible for n-nonane. If the %D for n-nonane is greater than 30, note the nonconformance in the case narrative. It should be noted that the %Ds are calculated when CFs are used for the initial calibration and percent drifts are calculated when calibration curves using linear regression are used for the initial calibration.

Actions:

If %RSD > 25% for target compounds or a correlation coefficient < 0.99, estimate positive results (J) and use professional judgment to qualify nondetects.

If % D > 25% (> 30 for nonane), estimate positive results (J) and nondetects (UJ).

CALIBRATIONS VERIFICATION

Compliance requirements for satisfactory instrument calibration are established to ensure that the instrument is capable of producing and maintaining acceptable quantitative data.

Date of initial calibration: _____ 06/22/16 _____

Dates of continuing calibration verification: _____ 06/30/15 _____

Dates of final calibration verification: _____ 06/30/16 _____

Instrument ID numbers: _____ GCDE _____

Matrix/Level: _____ AQUEOUS/MEDIUM _____

DATE	LAB FILE ID#	ANALYTE	CRITERIA OUT RFs, %RSD, %D, r	SAMPLES AFFECTED
Initial and continuing calibration meet method specific requirements				

A separate worksheet should be filled for each initial curve

DATA REVIEW WORKSHEETS

All criteria were met _____
 Criteria were not met and/or see below X

V A. BLANK ANALYSIS RESULTS (Sections 1 & 2)

The assessment of the blank analysis results is to determine the existence and magnitude of contamination problems. The criteria for evaluation of blanks apply only to blanks associated with the samples, including trip, equipment, and laboratory blanks. If problems with any blanks exist, all data associated with the case must be carefully evaluated to determine whether or not there is an inherent variability in the data for the case, or if the problem is an isolated occurrence not affecting other data. A Laboratory Method Blank must be run after samples suspected of being highly contaminated to determine if sample carryover has occurred.

List the contamination in the blanks below. High and low levels blanks must be treated separately.

Laboratory blanks

DATE ANALYZED	LAB ID	LEVEL/ MATRIX	COMPOUND	CONCENTRATION UNITS
---------------	--------	---------------	----------	---------------------

 METHOD BLANKS MEET THE METHOD SPECIFIC CRITERIA EXCEPT
 FOR THE CASES DESCRIBED IN THIS DOCUMENT.

 06/30/16 OP47988-MB Aqueous/low C11-C22_(Aromatics) 34.0 ug/L

Note: No action taken, blank concentration below the reporting limit. The laboratory qualified the results with a B qualifier.

Field/Trip/Equipment

DATE ANALYZED	LAB ID	LEVEL/ MATRIX	COMPOUND	CONCENTRATION UNITS
---------------	--------	---------------	----------	---------------------

 NO TRIP/FIELD/EQUIPMENT BLANKS ANALYZED ASSOCIATED WITH THIS
 DATA PACKAGE.

DATA REVIEW WORKSHEETS

All criteria were met X
Criteria were not met and/or see below

V B. BLANK ANALYSIS RESULTS (Section 3)

Blank Actions

The ALs for samples which have been diluted should be corrected for the sample dilution factor and/or % moisture, where applicable. Peaks must not be detected above the Reporting Limit within the retention time window of any analyte of interest. The hydrocarbon ranges must not be detected at a concentration greater than 10% of the most stringent MCP cleanup standard. Specific actions area as follows:

If the concentration is < sample quantitation limit (SQL) and < AL, report the compound as not detected (U) at the SQL.

If the concentration is \geq SQL but < AL, report the compound as not detected (U) at the reported concentration.

If the concentration is > AL, report the concentration unqualified.

DATA REVIEW WORKSHEETS

All criteria were met X
Criteria were not met and/or see below

SURROGATE SPIKE RECOVERIES

Laboratory performance of individual samples is established by evaluation of surrogate spike recoveries. All samples are spiked with surrogate compounds prior to sample analysis. The accuracy of the analysis is measured by the surrogate percent recovery. Since the effects of the sample matrix are frequently outside the control of the laboratory and may present relatively unique problems, the validation of data is frequently subjective and demands analytical experience and professional judgment.

List the percent recoveries (%Rs) which do not meet the criteria for surrogate recovery.

Matrix: solid/aqueous

SAMPLE ID	SURROGATE COMPOUND				ACTION
	S1	S2	S3	S4	

SURROGATE STANDARDS RECOVERIES WITHIN LABORATORY CONTROL
LIMITS.

S1 = o-Terphenyl 40-140%
S2 = 2-Fluorobiphenyl 40-140%
S3 = 1-Chlorooctadecane 40-140%
S4 = 2-Bromonaphthalene 40-140%

QC Limits (%)* (Aqueous)

LL to UL 40 to 140 40 to 140 40 to 140 40 to 140

QC Limits* (Solid)

LL to UL to to to to

Note: No action. % recoveries within laboratory control limits in second column.

It is recommended that surrogate standard recoveries be monitored and documented on a continuing basis. At a minimum, when surrogate recovery from a sample, blank, or QC sample is less than 40% or more than 140%, check calculations to locate possible errors, check the fortifying standard solution for degradation, and check changes in instrument performance.

If the cause cannot be determined, reanalyze the sample unless one of the following exceptions applies:

- (1) Obvious interference is present on the chromatogram (e.g., unresolved complex mixture);
- (2) The surrogate exhibits high recovery and associated target analytes or hydrocarbon ranges are not detected in sample.

If a sample with a surrogate recovery outside of the acceptable range is not reanalyzed based on any of these aforementioned exceptions, this information must be noted on the data report form and discussed in the Executive Report. Analysis of the sample on dilution may diminish matrix-related surrogate recovery problems. This approach can be used as long as the reporting limits to evaluate applicable MCP standards can still be achieved with the dilution. If not, reanalysis without dilution must be performed.

DATA REVIEW WORKSHEETS

All criteria were met _____
 Criteria were not met and/or see below __N/A__

VII. A MATRIX SPIKE/MATRIX SPIKE DUPLICATE (MS/MSD)

This data is generated to determine long term precision and accuracy in the analytical method for various matrices. This data alone cannot be used to evaluate the precision and accuracy of individual samples.

At the request of the data user, and in consideration of sample matrices and data quality objectives, matrix spikes and matrix duplicates may be analyzed with every batch of 20 samples or less per matrix.

- **Matrix duplicate** - Matrix duplicates are prepared by analyzing one sample in duplicate. The purpose of the matrix duplicates is to determine the homogeneity of the sample matrix as well as analytical precision. The RPD of detected results in the matrix duplicate samples must not exceed 50 when the results are greater than 5x the reporting limit.
- The desired spiking level is 50% of the highest calibration standard. However, the total concentration in the MS (including the MS and native concentration in the unspiked sample) should not exceed 75% of the highest calibration standard in order for a proper evaluation to be performed. The purpose of the matrix spike is to determine whether the sample matrix contributes bias to the analytical results. The corrected concentrations of each analyte within the matrix spiking solution must be within 40 - 140% of the true value. Lower recoveries of n-nonane are permissible but must be noted in the narrative if <30%.

MS/MSD Recoveries and Precision Criteria

Sample ID: _____ - _____ Matrix/Level: _____ - _____

List the %Rs, RPD of the compounds which do not meet the QC criteria.

MS OR MSD	COMPOUND	% R	RPD	QC LIMITS	ACTION

Note: No MS/MSD sample analyzed with this data package. Blank spike/blank spike duplicate used to assess accuracy. % recoveries and RPD within laboratory control limits. No action taken.

DATA REVIEW WORKSHEETS

All criteria were met X
 Criteria were not met and/or see below

No action is taken on MS/MSD results alone to qualify the entire case. However, used informed professional judgment, the data reviewer may use the MS/MSD results in conjunction with other QC criteria and determine the need for some qualification of the data. In those instances where it can be determined that the results of the MS/MSD affect only the sample spiked, the qualification should be limited to this sample alone. However, it may be determined through the MS/MSD results that the laboratory is having a systematic problem in the analysis of one or more analytes, which affects the associated samples.

2. MS/MSD – Unspiked Compounds

List the concentrations of the unspiked compounds and determine the % RSDs of these compounds in the unspiked sample, matrix spike, and matrix spike duplicate.

COMPOUND	CONCENTRATION		MSD	%RPD	ACTION
	SAMPLE	MS			

Criteria: None specified, use %RSD \leq 50 as professional judgment.

Actions:

If the % RSD > 50, qualify the results in the spiked sample as estimate (J).

If the % RSD is not calculable (NC) due to nondetect value in the sample, MS, and/or MSD, use professional judgment to qualify sample data.

A separate worksheet should be used for each MS/MSD pair.

DATA REVIEW WORKSHEETS

All criteria were met X
Criteria were not met and/or see below

VIII. LABORATORY CONTROL SAMPLE (LCS/LCSD) ANALYSIS

This data is generated to determine accuracy of the analytical method for various matrices.

1. LCS Recoveries Criteria

List the %R of compounds which do not meet the criteria

LCS ID	COMPOUND	% R	QC LIMIT	ACTION
--------	----------	-----	----------	--------

 LCS_RECOVERY_WITHIN_LABORATORY_CONTROL_LIMITS

Criteria:

- * Refer to QAPP for specific criteria.
- * The spike recovery must be between 40% and 140%. Lower recoveries of n-nonane are permissible. If the recovery of n-nonane is <30%, note the nonconformance in the executive narrative. RPD between LCS/LCSD must be < 25%.

Actions:

Actions on LCS recovery should be based on both the number of compounds that are outside the %R and RPD criteria and the magnitude of the exceedance of the criteria.

If the %R of the analyte is > UL, qualify all positive results (j) for the affected analyte in the associated samples and accept nondetects.

If the %R of the analyte is < LL, qualify all positive results (j) and reject (R) nondetects for the affected analyte in the associated samples.

If more than half the compounds in the LCS are not within the required recovery criteria, qualify all positive results as (J) and reject nondetects (R) for all target analyte(s) in the associated samples.

2. Frequency Criteria:

Where LCS analyzed at the required frequency and for each matrix (1 per 20 samples per matrix)? Yes or No.

If no, the data may be affected. Use professional judgment to determine the severity of the effect and qualify data accordingly. Discuss any actions below and list the samples affected. Discuss the actions below:

DATA REVIEW WORKSHEETS

All criteria were met _____
 Criteria were not met and/or see below N/A

IX. FIELD/LABORATORY DUPLICATE PRECISION

Sample IDs: _____ - _____

Matrix: _____ - _____

Field/laboratory duplicate samples may be taken and analyzed as an indication of overall precision. These analyses measure both field and lab precision; therefore, the results may have more variability than laboratory duplicates which measures only laboratory performance. It is also expected that soil duplicate results will have a greater variance than water matrices due to difficulties associated with collecting identical field duplicate samples.

COMPOUND	SQL	SAMPLE CONC.	DUPLICATE CONC.	RPD	ACTION
No field/laboratory duplicate analyzed with this data package. MS/MSD % recovery RPD used to assess precision. RPD within laboratory and validation guidance document criteria (+ 50 % RPD) for analytes concentration ≥ 5 SQL.					

Criteria:

The project QAPP should be reviewed for project-specific information.
 RPD $\pm 30\%$ for aqueous samples, RPD $\pm 50\%$ for solid samples if results are \geq SQL.
 If both samples and duplicate are < 5 SQL, the RPD criteria is doubled.

SQL = soil quantitation limit

Actions:

If both the sample and the duplicate results are nondetects (ND), the RPD is not calculable (NC). No action is needed.

Qualify as estimated positive results (J) and nondetects (UJ) for the compound that exceeded the above criteria.

If one sample result is not detected and the other is ≥ 5 x the SQL qualify (J/UJ).

Note: If SQLs for the sample and duplicate are significantly different, use professional judgment to determine if qualification is appropriate.

If one sample value is not detected and the other is < 5 x the SQL, use professional judgment to determine if qualification is appropriate.

DATA REVIEW WORKSHEETS

All criteria were met X
Criteria were not met and/or see below

XI. COMPOUND IDENTIFICATION

The compound identification evaluation is to verify that the laboratory correctly identified target analytes as well as tentatively identified compounds (TICs).

1. Verify that the target analytes were within the retention time windows.
 - Retention time windows must be re-established for each Target EPH Analyte each time a new GC column is installed, and must be verified and/or adjusted on a daily basis.
 - The n-nonane (n-C9) peak must be adequately resolved from the solvent front of the chromatographic run.
 - All surrogates must be adequately resolved from the Aliphatic Hydrocarbon and Aromatic Hydrocarbon standards.
 - For the purposes of this method, adequate resolution is assumed to be achieved if the height of the valley between two peaks is less than 25% of the average height of the two peaks.
 - The n-pentane (C5) and MtBE peaks must be adequately resolved from any solvent front that may be present on the FID and PID chromatograms, respectively.
- 1a. Aliphatic hydrocarbons range:
 - Determine the total area count for all peaks eluting 0.1 minutes before the retention time (Rt) for n-C9 and 0.01 minutes before the Rt for n-C19.
 - Determine the total area count for all peaks eluting 0.01 minutes before the Rt for n-C19 and 0.1 minutes after the Rt for n-C36.

Are the aliphatic hydrocarbons range properly determined?

Yes? or No?

Comments:

- 1b. Aromatic hydrocarbons range:
 - Determine the total area count for all peaks eluting 0.1 minutes before the retention time (Rt) for naphthalene and 0.1 minutes after the Rt for benzo(g,h,i)perylene.
 - Determine the peak area count for the sample surrogate (OTP) and fractionation surrogate(s). Subtract these values from the collective area count value.

Are the aliphatic hydrocarbons range properly determined?

Yes? or No?

Comments:

DATA REVIEW WORKSHEETS

All criteria were met X
Criteria were not met and/or see below

2. If target analytes and/or TICs were not correctly identified, request that the laboratory resubmit the corrected data.
3. Breakthrough determination - Each sample (field and QC sample) must be evaluated for potential breakthrough on a sample specific basis by evaluating the % recovery of the fractionation surrogate (2-bromonaphthalene) and on a batch basis by quantifying naphthalene and 2-methylnaphthalene in both the aliphatic and aromatic fractions of the LCS and LCSD. If either the concentration of naphthalene or 2-methylnaphthalene in the aliphatic fraction exceeds 5% of the total concentration for naphthalene or 2-methylnaphthalene in the LCS or LCSD, fractionation must be repeated on all archived batch extracts.

NOTE: The total concentration of naphthalene or 2-methylnaphthalene in the LCS/LCSD pair includes the summation of the concentration detected in the aliphatic fraction and the concentration detected in the aromatic fraction.

Comments: Concentration in the aliphatic fraction < 5% of the total
concentration for naphthalene and 2-methylnaphthalene

4. **Fractionation Check Standard** – A fractionation check solution is prepared containing 14 alkanes and 17 PAHs at a nominal concentration of 200 ng/μl of each constituent. The Fractionation Check Solution must be used to evaluate the fractionation efficiency of each new lot of silica gel/cartridges, and establish the optimum hexane volume required to efficiently elute aliphatic hydrocarbons while not allowing significant aromatic hydrocarbon breakthrough. For each analyte contained in the fractionation check solution, excluding n-nonane, the Percent Recovery must be between 40 and 140%. A 30% Recovery is acceptable for n-nonane.

Is a fractionation check standard analyzed?

Yes? or No?

Comments: Not applicable.

DATA REVIEW WORKSHEETS

All criteria were met X
Criteria were not met and/or see below

XII. QUANTITATION LIMITS AND SAMPLE RESULTS

The sample quantitation evaluation is to verify laboratory quantitation results.

In order to demonstrate the absence of aliphatic mass discrimination, the response ratio of C28 to C20 must be at least 0.85. If <0.85, this nonconformance must be noted in the laboratory case narrative.

The chromatograms of Continuing Calibration Standards for aromatics must be reviewed to ensure that there are no obvious signs of mass discrimination.

Is aliphatic mass discrimination observed in the sample? Yes? or No?

Is aromatic mass discrimination observed in the sample? Yes? or No?

1. In the space below, please show a minimum of one sample calculation:

MC46423-1 EPH (C11 – C22, Aromatics) RF = 124800

[] = (1965590)/(124800)

[] = 15.75 ppb Ok

MC46423-1 EPH (C19 – C36, Aliphatics) RF = 77820

[] = (1214579)/(77820)

[] = 15.61 ppb Ok

DATA REVIEW WORKSHEETS

2. If requested, verify that the results were above the laboratory method detection limit (MDLs).
3. If dilutions performed, were the SQLs elevated accordingly by the laboratory? List the affected samples and dilution factor in the table below.

SAMPLE ID	DILUTION FACTOR	REASON FOR DILUTION

If dilution was not performed, estimate results (J) for the affected compounds. List the affected samples/compounds:
