Microbial Commercial Activity Notice (MCAN) Submission

to the

U.S. Environmental Protection Agency Office of Pollution Prevention and Toxics Chemical Control Division New Chemicals Notice Management Branch

TS Number: J-012-1WK8NE

Date of Submission: March 21, 2012

Submitter:

Submitted to: TSCA Document Processing Center (7407)

Room L-100

Office of Pollution Prevention and Toxics


U.S. Environmental Protection Agency

1200 Pennsylvania Ave., N.W.

Washington, D.C. 20460

Certification of Information

I certify that to the best of my knowledge and belief: the company named in this submission intends to manufacture, import, or process for a commercial purpose, other than in small quantities solely for research and development, the microorganism identified in this submission. All information provided in the submission is complete and truthful as of the date of this submission. I am including with this submission all test data in my possession or control and a description of all other data known to or reasonably ascertainable by me as required by 40 C.F.R. § 725.160 or § 725.260.



SUBSTANTIATION OF CONFIDENTIALITY FOR THE INFORMATION CLAIMED AS CONFIDENTIAL BUSINESS INFORMATION IN THE MCAN

The following information is submitted in accordance with 40 C.F.R. § 725.94. A proposed generic name is provided for the microorganism identity in Section 1.4 of this MCAN pursuant to 40 C.F.R. §§ 725.80(a)(1) and 725.85(a)(3)(ii). A proposed category of use and a generic use description is provided in Section 1.5 of this MCAN pursuant to §§ 725.80(a)(2) and 725.88(b).

- A. The nature of the Company's business is relatively unique in U.S. and international commerce. The technology is such that a competitor would be able to discern the production of the microorganism if Submitter Identity, Microorganism Identity, Process Information, Use, and Internal Company Documents (excluding health and safety studies) are publicly disclosed. A competitor, upon discovering this information, would be able to manufacture, sell or use the microorganism with no investment in research or development, all to our company's detriment. Consequently, such a disclosure is intolerable.
- B. This information should be held in confidence indefinitely, until the technology is obsolete, or until the microorganism is widely known as a result of competing research.
- C. Our company has kept guarded information related to Submitter Identity, Microorganism Identity, Process Information, Use, and Internal Company Documents (excluding health and safety studies) so that others cannot discover the commercial utility of this information. Only those with a need to know have access to this information.
- D. Our company has disclosed information in these areas to outside legal counsel who appropriately protect its confidential status. Only those on a need-to-know basis are aware of the association between the organism and its location of manufacture and production so that its commercial utility cannot be discovered. No confidential information or licenses under existing patents and patent applications have been disclosed or granted to others without secrecy agreements.

- E. There are no advertising or promotional materials for the microorganisms. No Material Data Safety Sheet(s) is publicly available. No trade publications reference the microorganism. The Company's development of the microorganism has been held in strictest confidence such that no competitor is aware that this particular microorganism is in use. No Federal, local, or state agency or court has public files disclosing the Company's identity, process, or referenced internal documents in connection with the microorganism.
- F. The Company, pursuant to 40 C.F.R. § 725.92(b) and § 725.95(e), claims as confidential, any reference to the microorganism's identity as well as any information that would facilitate the discovery of its identity in: (1) health and safety studies conducted by the submitter; and (2) published scientific journal articles submitted with the MCAN. Disclosure of the microorganism's identity would reveal confidential processing manufacturing trade secret, unrelated to health and the environment. Less specific identity information is sufficient to interpret the references provided, because the results and conclusions of the researcher are fully disclosed by the articles.
- G. There are patents associated with the MCAN strain. However, the microorganism is only one of many microorganisms that have been disclosed in many of these patents. Thus, the identity of the microorganism, which is the subject of this MCAN, should be treated as confidential because these patents do not disclose the microorganism's identity, *per se*. Furthermore, the mere existence of such patents does not necessarily indicate that this microorganism or any other member of the category for which patent claims have been made, is utilized in U.S. commerce.
- H. No Federal agency or court has ruled on the confidentiality of the microorganism.
- I. A competitor, upon discovering this information, would be able to manufacture, sell or use the microorganism with no investment in research or development, all to our company's detriment. Consequently, such a disclosure is intolerable.
- J. EPA disclosure of the information claimed as confidential would allow a competitor to enter the market more easily, benefitting greatly from our research and development

expenditures. Competitors with the facilities, personnel and expertise to produce these microorganisms quickly, would benefit greatly at our expense, significantly reducing their own research and development time.

K.	The microorganism will not leave the site of production in a viable state thereby
	allowing for testing in a form that is accessible to the public or its competitors. The cost
	to a competitor, in time and money, to develop appropriate use conditions would be
	approximately The protection of
	the organism's confidentiality and secure handling will impede product analysis by
	others.

1.0 INTRODUCTION

1.1	Overview
	MCAN submission is for a specific genetically modified organism that will be used for
	Wier it submission is for a specific generically mounted organism that will be used for
	potential chemical uses.
The c	chemical substance produced
The c	chemical substance that is the subject of this request,
	Ongoing R&D efforts are conducted in compliance with the R&D exemptio
at 40	C.F.R §725.234 and 725.235

1




Figure 1a



Figure 1b.

1.2 Purpose of the MCAN



1.3 Contact Information

In accordance with 40 C.F.R. § 725.155(c), the following information is provided.





An agent letter is provided as **Attachment 1**.

1.4 Proposed Generic Name

The explicit and formal biological name of the microorganism
A generic name for this
microorganism that is in accord with 40 C.F.R. § 725.85 is "microorganism modified."
The manufacturer of the MCAN microorganism considers
as highly confidential the identity
This gene distinguishes the MCAN microorganism
from more conventional strains and contributes new and useful performance
properties to the microorganism. Nondisclosure of the specific gene used to modify the
microorganism is necessary to reduce the likelihood of a competitor manufacturing a similar
product without the considerable investment in research and development required to develop
such a product.

1.5 Proposed Use Category and Generic Use Description

The submitter proposes the following generic use description:


"chemical production." This description protects from disclosure the confidential process and purpose of the manufacture while disclosing, with respect to exposure and release, the chemical to be produced.

2.0 MICROORGANISM IDENTITY INFORMATION

In accordance with 40 C.F.R. § 725.155(d), the following information is provided.

2.1 Recipient Strain Identification –
The host organism,
does not specify biosafety levels for their microorganisms, they identify strains that have the
potential to be pathogenic. The parental strain is not identified as being potential pathogenic.
was taken through clonal (single cell) purification and given the Further information on the human health and environmental hazard assessment of is presented in
Sections 3.0 and 4.0 below.
and designated as Biosafety Level 1 (BSL-1).

2.1.1 Morphological differences



2.1.2 Genetic divergence

performed genotyping of	at the plastidic 23S rRNA
locus. As seen below in Figure 3, the genome sequences of these two s	species are clearly
different. The PCR method used to perform this genotyping can be use	ed to distinguish
species.	





Figure 3.





2.1.3 Biochemical/physiological differences


Figure 5

identical conditions were analyzed. The two strains exhibited very different growth rates, These differences are illustrated in Figure 5 and Table 1.	grown under
These differences are illustrated in Figure 5 and Table 1.	identical conditions were analyzed. The two strains exhibited very different growth rates,
	These differences are illustrated in Figure 5 and Table 1.

St	rain				
_					
Day					
8					
Day					
v.					
Day					



2.1.2.1 Relatedness of



Species per source collection	Date Acquired	Genotype (23S)	Source collection no.	Isolated from
			_	
			_	
			_	
			-	
			_	

Table 2.

The nucleotide sequences of the nine different 23S genotypes along with their accession numbers and species designation as per their respective source collection are provided in **Attachment 4** while their taxonomic relationship is diagrammed below in Figure 6.



Figure 6. Cladogram based on partial, plastidic 23S rDNA sequences

2.1.2.2 Genetic characterization

A complete genome sequence is available. We have carried out shotgun sequencing (Roche 454) and Illumina paired end reads as well as transcriptome analysis of the strain. The genome size of the organism is well as transcriptome analysis of the strain. The genome size of the organism is well as well as generate additional sequence information in order to resolve ambiguous assemblies and issues around highly homologous alleles and gene families through Roche paired end as well as BAC sequencing. We have also completed work aimed at characterizing the genome size of the organism through alternative methods including contour-clamped homogeneous electric field (CHEF) gel analysis and Fulgen staining (Table 3 and figure 7) followed by FACS (Flourescent Activated Cell Sorting) analysis and these results are in excellent agreement with our sequencing work. The organism appears to be a diploid.

Organism	Genome size by sequencing (Mbp)	Genome size by Fulgen staining (Mbp)

Table 3. Genome sizes as assessed by sequencing versus Fulgen staining.






Figure 7. Genome size




23S rDNA genotype



Figure 9. CHEF gel analysis of





All of the above data demonstrate that on multiple levels, there are no gross differences between

2.2 Features of the New Microorganism has an active strain improvement program to the wild-type strain To maximize the probability of generating mutants with the desired traits, both physical and chemical mutagens are used to exploit the differences in their modes of DNA alteration. Mutants that demonstrate increased productivities such are selected for subsequent rounds of strain improvement. Strain was taken through a classical strain improvement program utilizing both physical and chemical mutagens, and given the internal designation (Figure 8 above). Morphological assessment with brightfield and fluorescence staining has shown that strain (Attachment 5). It has also been determined that 2.3 Morphological and Physiological Features of New Microorganism We went on to transform the classically improved strain with a



Figure 11.

The resulting strain was given the	and is the
subject of this MCAN submission.	
has a similar morphology in terms of size and	compared to the
unmodified (Attachment 5) and now has the ability to	. Physiological
characterizations to determine similarities and any potential differences be	etween the original
parental strain	

2.3.1 Growth rates


As shown in Table 4 there are no significant growth rate differences in these strains.

Growth rates were determined by optical density measurements as detailed in **Attachment 6**.

Strain	Carbon Source	Growth Rate (h ⁻¹)



2.3.2 Growth temperatures



5






2.3.3	
As shown in Figure 13 below,	











2.3.5

Table 5.









2.4	Genetic Construction
2.4.1	The taxonomy of the donor organism
derive	is a stable transformant of strain, a classically improved strain d from The component parts of the transgene introduced into
	This is shown schematically in Figure 16 below.
	Figure 16
2.4.2	Description of traits for which the microorganism was selected

Additional information concerning the	
characteristics of the new strain can be found in Section 2.3 above.	
2.4.3 Genetic modifications	
We have developed the means to transform	
Transformation in is mediated by	
homologous recombination and as such, integrated cassettes have excellent stability.	
Construction of integrative vectors requires some knowledge of host genomic sequences	because
the vectors integrate through homologous recombination. Integrative vectors are design	ed such
that is cloned and modifie	d if
necessary to contain a suitable restriction site at the approximate midpoint of the genome	ic DNA.
Heterologous genes of interest can then be inserted into this site and the entire cassette the	hen
cloned into a suitable vector for replication in E. coli, such as pUC19 or pBlueScript. The	he
resulting cassette consists of the heterologous gene now flanked on either side by what v	was
formerly a contiguous piece of host genomic DNA, cloned into an E. coli-based vector.	
Prior to transformation, the DNA is linearized with restriction enzymes such that the	
heterologous gene of interest, flanked on either side by genomic DNA, is physical	ılly
separated from the rest of the plasmid sequences necessary for maintenance in E. coli.	
Transformation of the linearized construct allows for the insertion of the heterologous ge	ene at the
precise location in the genome defined by the flanking sequences through homological	ogous
recombination (Figure 16 above).	
Codon usage and gene optimization in	
Proper expression of transgenes requires that the codon usage of the transgene matches t	the
specific codon bias of the organism in which it is being expressed. The codon usage in t	the
recipient strain was deduced from extensive analysis cDNA sequences. A	All genes
encoding proteins in this document have been codon optimized to ensure efficient translation	ation of
the resulting mRNA. A sequence comparison of the wild type sequence compared	to the
codon optimized sequence is provided in Attachment 11 .	

2.4.5 Detailed description of the genetic construction


Figure 17 below is the complete nucleotide sequence of the transforming DNA, including its		
vector backbone, used to generate		

Figure 17. Complete coding sequence	
Γhe	
Flanking sequences are derived from the	

2.4.6 Stability of Inserted Genetic Material				
All modifications were integrated directly into chromosomal DNA. No loss of any inserted				
genetic material has been observed. The				
was demonstrated to be 100% stable as follows. The primary transformant giving rise to				
generated as described earlier and				
was grown in liquid culture with as the sole carbon source (non-selective conditions).				
After 30 generations, cells were plated to solid medium containing as the sole carbon				
source. 48 individual colonies were then picked from this plate and transferred to a 96 well plate				
made up with medium containing as the sole carbon source. Forty eight of forty eight				
clones transferred exhibited growth on the				
generations of growth in the absence of selection, indicating the inserted transgene is				
100% stable. To further validate the stability of the transgene and its correct targeting to				

Figure 18.





.4.7 PCR Amplification
figure 20 shows that the primer pairs are only capable of amplifying a product in the transgenic
ne and it is of the correct size. The resulting PCR product was then cloned and subjected to
equencing, the results of which are shown in Figure 21.

Figure 20. PCR analysis



3.0 ECOLOGICAL ASSESSMENT

3.1 Phenotype

The strain phenotype request is understood to refer to the expression of the genes of the
organisms as well as the influence of environmental factors and random variation. The
interaction between these factors may be represented as genotype + environment + random
variation \rightarrow phenotype. The strain phenotype produces
3.2 Habitat, Geological Distribution and Source of the Recipient Microorganism
3.2.1 Description of
Regarding the origins





3.2.2 Geographical distribution

3.2.2.1 Overview



	İ
	i
3.2.2.2 Aquatic environments	
3.2.2.3 Sewage environments	
	1



3.2.2.4 Food contact



Table 6.

Although found in food products and thereby ingested, laboratory experiments in rodents and primates have shown that are unharmed by the digestive process and pass
through the gut into the feces with little or no reproduction or harm to the host
3.3 Detection, Survival, and Dissemination under Relevant Environmental Conditions
3.3.1 Control technologies, inactivation protocols, and fermentation conditions
is grown in a contained environment in a production facility, and it is highly unlikely for
any viable to be present in the final product, as
inactivate all viable cells. The organism is classified as a BSL-1 and Good Large
Scale Practice (GLSP) organism (GLSP at >10 L). Control technologies at release points and
inactivation protocols are in place to reduce or eliminate the dissemination of viable modified
microorganisms. Please refer to Section 7.0 for further details on the air, solid waste, and
wastewater inactivation treatments and containment descriptions. In addition to inactivation
protocols and control technologies, as discussed in Section 7.3
the fermentation process engineering controls and protocols, restrict the ability to
reproduce, disseminate, and spread outside the contained vessel.
reproduce, disseminate, and spread edistice the contained vesser.

3.3.2 Description of method for detecting the microorganism in the environment

Section 2.4.7 of the MCAN provides details concerning a method for detecting the microorganism in the environment utilizing PCR analysis.

3.4 Anticipated Biological Interactions with Target Organisms and Other Organisms

	Host range: The range of host species or cell types which the modified production strain is
	able to infect or parasitize is expected to be no different than that for the recipient strain that
	is typically used in chemical production today
•	Target organism: The use of the production strain is
	is not designed to act upon a particular organism during the production process or otherwise.
•	Competitors: During fermentation, competing microorganisms include normal airborne
	microorganisms such as that are considered to be contaminants. If
	left unchecked, they will compete with the production organism for nutrients and resources,
	which will cause stress to the modified microorganism, which will reduce yields
	Therefore, aseptic techniques are employed during the whole process from seed
	propagation to fermentation.
•	Prey: The host and modified microorganism is not designed to prey upon living organisms
	as a food source.
•	Hosts: The use of the production strain is for
	microorganism is not designed to be a host or to infect or feed upon another living organism.
•	Symbionts: The production strain is not designed to be an organism in a symbiotic
	relationship. It is not designed to serve as a host in which the presence of a smaller symbiont
	beneficiary would be present.
•	Parasites: No significant interactions with parasites were discussed in a review of publically

Pathogens: No significant interactions with pathogens were discussed in publically available

scientific literature. The production strain or host strains are not designed to enhance any

pathogens such as Escherichia coli O157, Clostridium botulinum, etc.

available scientific literature.

3.4.1 Effects on plants

There are no published studies on the effects of the recipient microorganism of this MCAN, on plants. However, there are studies available that describe the effects of that we have used as surrogate data to evaluate plant pathogenicity.
Published evidence indicates that several different species of plants serve as a natural habitat for with no reports of adverse effects. Therefore, species in the genus are unlikely to be a plant pathogen.
3.4.2 Effects on animals
Although the presence of

As noted, several animal studies demonstrate
In addition, we have summarized other animal studies in the human health hazard assessment section 4.0.
3.4.3 Potential for gene transfer
We were unable to locate any published articles discussing the potential for gene transfer
between and other organisms. Furthermore, we were unable to locate any published
articles on laboratory studies related to gene transfer between and other organisms.


Therefore, we have based our analysis on potential for gene transfer between microorganisms in
general. Microorganisms in the wild are constantly being exposed to DNA with an estimated
0.2-19 μg DNA/L in marine aquatic environments, 0.5-25.6 μg DNA/L for freshwater
environments, 1 mg/DNA/g of freshwater sediment, and 80-85 µg DNA/g of soil
DNA, although abundant in the environment, is not stable due to
chemicals and nucleases present in the environment, and DNA may also bind to other solids.
The half-life of DNA in wastewater is 1-13.8 minutes, 4.2-5.5 hours in freshwater, 3.4-83 hours
in marine waters, and 9.1-235 hours in sediments and soils
In addition, the recipient microorganism needs to be in a competent state and under ideal
environment conditions (i.e., temperature shifts, specific cautions, or electrical pulses) in order
for gene uptake to occur, and most prokaryotes, which are typically the group of microorganisms
to participate in genetic transfer from the environment, are not competent in the wild
Therefore, it is not surprising to note that despite the
detection of transgene from genetically modified plants in the soil, no gene transfer has been
documented to occur to any of the soil microbial population
Uptake of genes from the modified requires the host (recipient)
microorganism to have proper intracellular DNA stabilization mechanisms, accurate expression
of the gene, and appropriate posttranslational modification mechanism
Furthermore, the intergeneric gene inserted into the modified
3.4.4 Interactions in the environment – e.g., biogeochemical cycling, trophic level,
competitiveness, non-targets, etc.
3.4.4.1 Recipient microorganism
There is very little information in the public literature that specifically addresses the ecological
role of . role in the environment is as
. As

3.4.4.2 MCAN microorganism

The modified strain expresses a	
Based on Table 4, the growth rate of	
Therefore, growth characterizations of the modified strain show that is unlikely to grow faster than the unmodified strain.	it
(Figure 15).	

4.0 POTENTIAL HUMAN HEALTH HAZARD ASSESSMENT

has conducted a thorough assessment on the safety of both the host organism and
the genetic elements used to generate
are classified as Biosafety Level 1 (BSL-1) and have no
known hazard to laboratory personnel. These strains are non-pathogenic, harmless organisms.
Standard microbiological practices consistent with BSL-1 organisms for production volumes at
or below 10 L should be followed when working with this strain (Biosafety in Microbiological
and Biomedical Laboratories, U.S. Department of Health and Human Services, Centers for
Disease Control and Prevention, and National Institutes of Health, February 2007). For volumes
greater than 10 L, GLSP level treatment is appropriate for this organism. Both the facilities in
have levels of release controls beyond the requirements
of GLSP such as
Based on available scientific and patent literature, are safe. The
genetic elements used to generate are well characterized, come from well-known and
studied microorganisms, and have no known risks.
In addition, since
in dedition, since
·
has performed a significant number of stirred tank fermentations and
shake flasks experiments using
an Institutional Biosafety Committee (IBC) comprised of technically
qualified individuals, which has not observed or reported any adverse health effects during the
production of



with the host strain despite its practice of doing so if it has such information.

Extensive small–scale and large-scale cultivation of
has not generated any negative health or
environmental effect data. is not aware of any information on health effects that
accompany any EPA rule or order issued under TSCA section 4, 5, or 6 that applies to
Based on this assessment, concluded that cultivation of has no known
health risks.
4.1 History of Safe Use
4.1.1 Recipient microorganism
This species,
has not yielded any known or observed negative
effects on the safety or health of personnel.

4.1.3	Donor	organism
T.1.J	Donor	or gamem

4.2 Potential for human pathogenicity
The recipient strain has been classified as a BSL-1 organism by ATCC which means that it is not known to cause disease in healthy adult humans. Our review of primary research articles below concurs with the BSL-1 assessment that is a rare and opportunistic pathogen and not a vector for pathogens.

4.2.3 Pathogenicity/toxicity

_

4.2.4 Submission microorganism – effect of genetic monopathogenicity potential, behavior, effects of productions of the pathogenicity potential, behavior, effects of genetic monopathogenicity potential, behavior, effects of genetic monopathogenic monopath	
The modified	
compared to the unmodified recipient strain	This is not anticipated to
ncrease the pathogenicity potential of the modified	The growth rate of the modified
strain is the same or less than both the original parer	
4.3 Allergenicity	
There are no nublished remorts of any mambers of the	as a potential allergen.
There are no published reports of any members of the However, as is taxonomically related to	as a potential allergen.
and there are published studies on the allergenic potential of	of we have based our
	for consideration. In a study that
	k tests and a conjunctival provocation
rest concluded a higher incidence of a positive reaction to	could be found in children
who were already sensitive to mold or other allergens	
·	vas a "weak allergen" that may be
clinically important to populations sensitive to other allerg	•

4.4 Virulence

The ability for an organism to be pathogenic is promoted by the presence of virulence factors, such as proteases, lipases, and toxins. Because the gene modification to the organism was not shown through a literature search to be toxic or yield different toxicological results from the parental strain, surrogate information on the recipient strain is offered for the purpose of evaluating the anticipated behavior of the production strain.

4.5 Immunologic Reactions

The inserted genetic elements in this case do not appear to possess any intrinsic hazard potential given the history of safe use of the enzyme, and the gene modification to the organism was not shown through a literature search to be toxic or yield different immunological or toxicological results from the parental strain.

4.6 Antibiotic Resistance (ABR)

The final production strain does not contain any ABR genes, and therefore the modified
microorganism is as susceptible to antibiotics and to anti-fungals as the unmodified
microorganism
A PubMed and a Google Scholar search did not produce any articles related to resistance to
antibiotics when we conducted a search using the terms in
conjunction with metal, pesticide, or herbicide. A null result was expected since
is not known to have any effect upon the survivability of the

Since the inserted genetic elements in this case do not appear to possess any intrinsic ability to confer antibiotic resistance, data are being provided for the species in general based on the rationale that the gene modification to the organism has not been shown, through a literature search, to impact antibiotic resistance *per se*.

intergeneric sequence, it is reasonable to conclude that the modified strain is not expected to be
any different from other well-known
strains commonly found in nature.
4.7 Action as a Vector for Pathogens
A PubMed and Google Scholar search using the terms
pathogen* and "vector" does not return any articles, demonstrating that the presence of the
alone is unlikely to permit the production strain to act as a vector of pathogens.
A null result was expected since
There are
no studies that the submitter could locate that would indicate the donor strain itself,
acts as a vector for pathogens. Because the gene modification to the organism was not shown
through a literature search to yield different results from the parental strain, surrogate
information on the recipient strain is offered for the purpose of evaluating the anticipated
behavior of the production strain. The production strain is not expected to act as a vector of any
pathogen such as Escherichia coli or Clostridium botulinum.

Based on the absence of demonstrated adverse effects for the parental strain and for the inserted

- 5.0 MANUFACTURING PROCESS, USE, AND DISPOSAL OF THE STRAIN
- 5.1 Maximum Amount to be Produced in Year One and Year Three (by volume) and Estimation of Viability (cells per unit volume or cfu)

The modified strain will be used to produce in the United States. The expected volumes are as shown below in Table 7.

YEAR	# of FERMENTORS	MAX # BATCHES	UTIL RATE	# of BATCH/ Year	CFU per mL BROTH	CFU per BATCH	CFU per YEAR (CFU)	Total fermentor VOL per BATCH (L)	Average fermentation time/batch hours

Table 7. Estimated production volumes for a three-year period. Please refer to Attachment 12 for the laboratory protocol.

Utilization rate ("UTIL RATE") refers to the percentage of time during that year that the
Submitter anticipates that the manufacturing facilities will be operational. To capture the most
accurate total viable cell counts,
5.2 Manufacturing Process Description
5.2.1 Operational description

The are grown in closed fermentation and processing plant for are grown in closed fermentors, and the entire site provides containment for the use of such microorganisms through the use of berms, trenches and sumps which channel all liquid to

tanks where the cells can be inactivated by heat or chemical means.


Microorganisms are handled at the appropriate level of biosafety (BSL-1, GLSP) as designated by CDC/NIH guidelines. Handling of effluent streams, solid, liquid, and gaseous, are all designed to provide the appropriate reduction and control of the microorganism to minimize any off site releases.

The identity of the site at which the operation will occur is
, for the fermentation and subsequent post-fermentation
processing
5.2.2 Transportation
Master cultures are stored in
Maintenance and storage of the master culture stock will be at the Submitter's
laboratory under the supervision of a qualified laboratory microbiologist. Cell propagation,
growth,
To transport the master culture stock
vials will be placed in a double zip-lock bag within a Styrofoam container. The
container will be filled with dry ice and transported via commercial shipper (e.g. Federal Expres
or equivalent) in accordance with Department of Transportation Hazardous Materials
requirements. The container will include a material safety data sheet (MSDS) that describes the
nature of the contents (Attachment 29). Instructions given on the MSDS describe containment
methods for an accidental release and disposal of the material. Upon receiving a master stock
culture, a staff microbiology technician at the facility, using standard laboratory procedures, wil
propagate the culture and prepare a working stock culture.

5.2.3 Cell propagation



Figure 24. Cell propagation diagram.



	The connection is then sterilized by
steam after the flask hose is attached.	
will be cleaned through chemical inactivation	methods such as bleach, and
autoclaved for the next use. Bleach has been shown to be effe	
autociaved for the next use. Breach has been shown to be circ	serive in macrivation

5.2.4 Fermentation

All process water from cleaning and related activities will be collected and sent to an on-site waste treatment area, where cells will be treated using heat, chemicals, or other suitable means to ensure inactivation.

All samples taken during the fermentation stages will be taken from a sample port and collected in sterile containers. The containers will be taken directly to the facility's laboratory for analysis. Broth will be purged from the ports prior to sampling and sent to an on-site waste treatment area, which will inactivate any residual cells through heat or chemical treatments.



Figure	


5.2.4	I		

released into the atmosphere by dryer exhaust fans.	. The remaining air will be
5.2.5	
5.2.6	

5.2.7 Waste disposal and sterilization procedures

will be
sterilized in the autoclave for re-use. With respect to sterilization of the transfer hose or lines
connecting the fermentation tanks, these will also be steam-sterilized.
After a tank has been emptied, the tank must be rinsed with water at the conclusion of the run.
The rinse water will then be sent to an on-site waste treatment area, where cells will be treated
using heat, chemicals, or other suitable means to ensure inactivation. After rinsing the tank, it
will be heat-sterilized with steam The tank sterilization will be
completed before any subsequent runs.
with liquids
being sent to an on-site waste treatment area.
the transfer hoses will all be flushed with water after use. The
inse water will then be sent to an on-site waste treatment area.

5.3 Product Description






Table 8:



6.0 WORKER EXPOSURE INFORMATION

The manufacturing site will be under the control of the Submitter. The production strain is
currently intended for use in the
which will meet the containment criteria to qualify for the full Tier 1 exemption. Future
production may expand to include other facilities that will also meet the containment criteria to
qualify for the full Tier 1 exemption. Workers will be grouped into shifts consisting of
people Lab Technician, Production Technicians, Operator/Mechanic, and Shift
Supervisor (Table 9 below). Shift teams provide coverage 24 hours/day, 7 days/week.
Additional personnel may be on shift, depending on workload.
Lab Technicians are responsible for the storage of working stock cultures, preparation of
inoculum (flask growth stages), and analyses of fermentation samples. Production Technicians
are responsible for operation of all manufacturing steps,
Maintenance Operator/Mechanics are
primarily responsible for addressing maintenance issues which arise during operations,
coordinating with Plant Site Maintenance Department for assistance as needed, and assisting the
production operators with plant operations. Engineers/Tech Managers and Roving Technicians
provide technical resources and supervision of all activities that occur on their shift

Position	Activity	PPE/Eng Control	Max # workers involved in each activity	Max duration Hours/day Days/year
Production Technicians	Sampling, Cleaning fermentors Inoculating	Gloves, safety glasses, lab coats		Actual exposure limited by engineering controls, PPE.
"Roving Techs"	Process support		_	Handling of cultures while sampling less than 1 hour/day
Lab	Culture	Gloves, safety		8 hour shifts
Technicians	transfer, shake flask	glasses, lab		250 days/year
	inoculation	Cours		Handling of cultures while sampling less than 1 hour/day
Engineer/Tech	Process	Gloves, safety		8 hour shifts
Manager	oversight	glasses, lab coats		250 days/year
Maintenance	Repair of	Gloves, safety		8 hour shifts
Operator	fermentation	glasses		250 days/year
Mechanics	equipment:			
	valves, pumps,			
	etc.			

Table 9: Worker exposure.

Exposure will be negligible, as the fermentation process is conducted in a totally closed system from which samples are drawn approximately once every 2-4 hours for in-process analysis. This sampling will only take about 1 minute. The Lab Technician is in charge of in-process sampling, but Production Technicians may also take samples.

In addition, facility operations are tailored to minimize exposure of personnel to viable microorganisms and to restrict releases of viable microorganisms into the environment. Employees will all be trained in aseptic techniques and in the potential health risks of microorganisms. Further, access to the fermentation lab and production areas – indeed, to the entire plant site – is strictly controlled and limited to employees and authorized visitors. Microorganisms are only handled in a manner intended to safeguard personnel health, through standard PPE and procedures.

7.0 POTENTIAL RELEASE INTO THE ENVIRONMENT

Cells from air exhaust, liquid, and solid waste are captured through extensive control technologies and inactivation protocols as described in Table 10 below.

Release	New substance released CFU/Day to Environment/ Control Technology	Media of Release	Control Technology	Efficiency
Process Wastewater	The site will inactivate viable cells, collect and treat wastewater, Minimal/no release to the environment.	Waste rinse water from cleaning of fermentors, transfer lines.	Heat or chemical inactivation in treatment tanks.	6 log reduction or greater, based on samples taken after treatment.

Solid waste	Daily plates,	Agar plates and	Autoclaved	6 log
from lab	and other solid	other solid lab	to ensure	reduction or
	waste.	waste	sterility.	greater based
		contaminated with the culture.		on heat inactivation
	no release to	the culture.		tests of
	the			biomass after
	environment.			drying.
	All solid waste			
	from labs			
	autoclaved.			
Exhaust from	Minimal/no	Air exhaust from		
fermentors	release to the	fermentors		
	environment.			
	liquid waste,			
	which is treated before release			
	before felease		capture and	
			condensate	
			out potential	
			aerosols into	
			wastewater,	
			where they will be	
			treated with	
			heat or	
			chemicals	
			before	
			release.	
Lab sink waste		Sample bottles	All samples	6 log

day 5 days /week from in process sampling. no release to the environment. All waste from labs autoclaved or chemical treated.	washed in lab sink going to drains that do not connect to treatment tanks.	are autoclaved or treated with chemicals to ensure inactivation before being washed in sinks.	reduction or greater based on chemical and/or heat inactivation data.

Table 10. Releases and control technology. Cells from air exhaust, liquid, and solid waste are captured through extensive control technologies and inactivation protocols.

7.1 Air Release

7.1.1 Fermentation

Air exhaust from the
vapor droplets (see Figure 24 Cell Propagation Diagram). Steps have been
taken to minimize this exposure.
a vent
allows air and smaller vapor particles to escape, while being designed to restrict liquid releases
via the exhaust air from the structure. The liquid releases are captured in a liquid stream, which is sent to the on-site wastewater treatment plant. With these controls in place, expected emissions
are estimated to be insignificant.
are estimated to be insignificant.
7.1.2

7.2 Wastewater Release

will produce lab scale wastewater from washing vials and flasks that come into contact with Lab ware containing viable cells is chemically treated or autoclaved before washing to ensure inactivation.
Process wastewater is generated in the fermentors and tanks collecting liquid vapor from The wastewater produced will be handled on-site at a wastewater treatment plant (WWTP) in accordance with a wastewater discharge permit from the The facility will discharge according to the permit requirements
Treatment steps included in the WWTP serve to restrict the release of viable into the environment. Critical elements of the wastewater treatment system that support this determination include:
 Heat or chemical treatment to inactivate organisms; Neutralization of the wastewater to pH level 6.2-9.0; Removal of floatable and settable solids through flocculation and secondary clarifiers. A portion of this sludge is then returned to the WWTP's aeration basins; another portion is sent back to the wet well for reprocessing; and the remaining 8-15% is collected, dewatered and transported off-site for incineration; and Discharge of clarified effluent after combining with storm water and non-contact coolant water in accordance with the WWTP's permit conditions.
Because of these treatment steps, and our knowledge that the modified at issue is neither will be restricted in effluent discharge.

7.2.1 Failed fermentation broth
Wighla will exist in the formantation heath. Thus, in the arrest that the autim contents of a
Viable will exist in the fermentation broth. Thus, in the event that the entire contents of a
fermentation vessel must be discarded due to a "failed" run, the contents will be inactivated by
heating to greater prior to discharging the broth to the WWTP. The
heat step results in at least a 6 log reduction in total viable cells as our data demonstrates greater
than 7 log reduction when cells are treated at a lower temperature
(Attachment 31).
7.3 Solid Waste
Solid waste from lab activities will be inactivated with chemicals or heat (autoclaving) to ensure
viable cells are not present in lab waste. Please refer to Section 5.2.7 entitled "Waste Disposal
and Sterilization Procedures" for details regarding the disposal of inactivated solid waste.
7.3.1 drying
Viable will not survive the on-site drying process. The final fermentation broth containing
the viable will be exposed to temperatures ranging from
, conditions that will inactivate the cells (see Attachment
32). Similarly, final fermentation broth to be sent off-site for drying will also go through a heat
treatment step: the fermentation broth will be heated to greater minutes, which
are conditions that will inactivate the cells (Attachment 31). These protocols are designed so
that viable will not survive the drying process, and therefore viable will not be
present in the
is then sent to landfill in accordance to local regulations and authorities.

7.4 Dissemination and Spread of Viable

The fermentation process engineering controls restrict the functional ability of
reproduce, disseminate, and spread. Specifically, the conditions in the fermentor are designed to
hold the organism at a physiological state at which they
The entire fermentation takes place an enclosed vessel, with several stages of containment to
prevent offsite releases. The threat of viable spreading into the environment during the
fermentation process is negligible.

As discussed in Section 7 above, cells from air exhaust, liquid, and solid waste are captured through extensive control technologies and inactivation protocols.

8.0 PROCEDURES FOR DISPOSAL OF ARTICLES, WASTE, CLOTHING, AND EQUIPMENT

All laboratory waste known to have been in contact with active biological materials are disposed of in designated biohazard waste receptacles. Biohazard waste is collected and killed using an autoclave heated to >121° C. Once deactivated this waste is disposed of via a private garbage handler with normal waste streams.

Lab clothing is collected and cleaned by an industrial service. Clothing also can be cleaned using a normal household detergent.

All other equipment, including bench tops are disinfected using commercial antimicrobial cleaner or equivalent after coming into contact with a modified microorganism. When acceptable, equipment will be placed into an autoclave and heated to >121° C in order to decontaminate prior to re-use.

9.0 SPILLS AND EMERGENCY PREPAREDNESS MEASURES

Good Microbiological Lab Practices and Good Large Scale Practices for BSL-1 and GLSP will be employed. Personnel routinely utilize the following practices in their normal operation to minimize exposure of personnel to viable microorganisms, and to restrict release of viable microorganisms to the environment. All of these practices must be utilized while conducting processes involving transgenic strains:

- 1. Employees are trained in aseptic techniques and in the potential health risks of microorganisms utilized so that potential biohazards can be understood and appreciated;
- 2. Access to fermentation lab and production areas, and the entire Plant site, is strictly controlled and limited to employees and authorized visitors;
- 3. Strains are handled in facilities intended to safeguard personnel health during handling;
- 4. Work surfaces where strains are handled are cleaned and decontaminated at least once per day with one or more of isopropyl alcohol, ethanol, Vesphene IISE, LpH se, or equivalent;
- 5. Contaminated wastes are either disposed via the lab sink which drains to the holding tanks (*e.g.*, broth samples after testing has been completed), or autoclaved or bleached prior to disposal (*e.g.*, used pipettes, used sample cups, used agar plates, empty cryovials);
- 6. Mechanical pipetting devices are used; mouth pipetting is prohibited;
- 7. Eating, drinking, smoking, and applying cosmetics are not permitted in the work areas. Food may only be stored in cabinets or refrigerators designated and used for this purpose only;
- 8. All persons wash their hands after they handle materials involving live strains;
- 9. In the interest of good personal hygiene, facilities (*e.g.*, hand washing sink, shower and change room) and protective clothing (*e.g.*, uniforms, lab coats, and safety glasses) are provided;

- 10. An insect and rodent control program is in effect;
- 11. The lab has windows that cannot be opened; the lab is air-conditioned;
- 12. Sample collection, transfer of culture fluids within/between systems, and processing of culture fluids are conducted in a manner that minimizes release of aerosols and maintains employee exposure to viable cells at a level that does not adversely affect the health and safety of employees;
- 13. A closed system which contains, or has contained, viable cells (e.g., a fermentation vessel) shall not be opened for maintenance until it has been thoroughly washed and sterilized by heat, chemicals or equivalent;
- 14. The site has institutional codes of practice that have been formulated and implemented among all employees to assure adequate control of health and safety matters; these include an industrial hygiene program, safety program and procedures, training, and periodic refresher training; and
- 15. The site has an emergency response plan, which includes provisions for handling spills.

10.0 PROCEDURES FOR TERMINATING THE ORGANISM

All fermentation broth is intended for batch may perform poorly or become accidentally contaminated with a competing organism. Cross-contamination incidents will be unlikely due to SOPs, including cleaning of the fermentors before and after use. However, in the event of cross-contamination or for other contamination reasons, termination may be desired. In such a case, the fermentation broth and all other liquid in the process would be processed with heat inactivation at temperature and time treatments or chemical treatment known to inactive the modified microorganisms. The fermentation equipment would undergo CIP rinses, heat sterilization and all wastewater and CIP liquids are treated to inactive the modified microorganism. Inactivation can be confirmed by samples plated onto growth media to show no viable are present in the treated wastewater.

11.0 HEALTH AND SAFETY DATA

As required by 40 C.F.R. §725.160, all testing regarding the health and environmental effects conducted on the microorganism known to the Submitter have been provided with this MCAN. Copies of the studies and references in the Submitter's possession and control have also been provided with this MCAN.

12.0 PROCESS DECRIPTION OF SITES NOT CONTROLLED BY THE SUBMITTER

The production strain is currently intended for use
which will meet the containment criteria to qualify for the
full Tier 1 exemption. However, future production may expand to include other facilities that
will also meet the containment criteria to qualify for the full Tier 1 exemption.