



2425 New Holland Pike, PO Box 12425, Lancaster, PA 17605-2425 • 717-656-2300 Fax: 717-656-2681 • www.lancasterlabs.com

Type I Data Package

Prepared for:

Olin Corporation Suite 200 3855 North Ocoee Street Cleveland TN 37312 CHECKED FOR COMPLETENESS OF PARAMETERS, ORDERED BY:

Clurecas

Project: Olin Wilmington, MA Superfund Site/6107090016 Soil Samples Collected on 06/08/11

SDG# OLN72

GROUP 1250628 SAMPLE NUMBERS

6310728-6310731

PA Cert. # 36-00037 NY Cert. # 10670 NJ Cert. # PA011 NC Cert. # 521

TX Cert. # T104704194-08A-TX

Through our technical processes and second person review of data, we have established that our data/deliverables are in compliance with the methods and project requirements unless otherwise noted or previously resolved with the client.

Authorized by:

Koma on Kayfman.

Date: 07/06/2011

Dana M. Kauffman Manager

Any questions or concerns you might have regarding this data package should be directed to your client representative, Nicole Maljovec at Ext. 1537.

Total Number of Pages 6



Table of Contents for SDG# OLN72

1.	Sample Reference List 1
2.	Analysis Request, Field Chain-of-Custody Record 2
3.	Methodology Summary/Reference 5
4.	Analysis Reports 6
5.	Hydrazines by LC/MS/MS Data
	a. Case Narrative-Conformance/Nonconform. Summary 14
	b. QC Summary 17
	c. Sample Data 22
	d. Standards Data 29
	e. Raw QC Data 51
	f. Preparation Logs 57
6.	Moisture Data



Sample Reference List for SDG Number OLN72 with a Data Package Type of I 12670 - Olin Corporation

Project: Olin Wilmington, MA Superfund Site/6107090016

Lab	Lab	
Sample	Sample	
<u>Number</u>	<u>Code</u>	Client Sample Description
6310728	448-D	OC-SS-448-0.0/1.0-DUP Grab Soil
6310729	448-1	QC-SS-448-0.0/1.0-XXX Grab Soil
6310730	448-1	OC-SS-448-0.0/1.0-XMS Grab Soil
6310731	448-1	OC-SS-448-0.0/1.0-MSD Grab Soil

15-8210180/820021/01261A

Client: Olin (Address: 3855	Olin Corporation 3855 North Ocoee St. Suite 200 Cleveland, TN 37312	te 200	Client	Client Project #: Work Site ID: Reports Sent To	#: 510 Will To: Ster	Client Project #: 5107090016 Work Site ID: Wilmington, MA Reports Sent To: Steve Morrow	∀ 3	1		Company Na	NVOICE INFO Company Name: Company Contact	e: Olin Corp	Corp	NVOICE INFO Company Name: Olin Corp Company Contact: ERG Accounts Payable	Shaded Areas for office use only	o esn eouro	ylk	
Phone: 423.: Raquested Turns Standard	Phone: 423-336-4511 Fax: Requested Turnersund Time (Shecility) Standard Rush (Lab Approval	A23-335-1466 Email SGMorrow@ Regulatory Programs: MADEP MCP Sul Report Requirements Level IV Package EDD Requirements: MACTEC EQUIS EZ EDD	5-1466 Program teiremen	Email rs: MADE rts Level N	<u> </u>	GMorrov			Email Rpt:	Address: Phone: Job #	** **	Sam	Same as Cllent Er	Cllent Email Quote# PO#				
M	-		(D) or Grab (G)	· · · · · · · · · · · · · · · · · · ·	033 - GC/NPD) c · ·					1 .	C/WG/WG)	NMH, UDMH ≥ ⅓		< Proservative Type (4) <-Bottle Type (5) Bottle Type (6)	Lancestr Lab	to La	91	
Sample ID OC-SS-448-0.0/1.0-DUP	DateTime Collected 1.0-DUP 5/35:00 AM	Fraction (1)	Sample Matr Composite (i	A0906) 3+13 (Cr+6 (7199)	DMF (Mod &	Perchorate (A janizarbyH → 3168 boM)			Ship	Comments (Special Instructions)	In Instruction	<u> </u>	7
						S	pecial	Instru	ctions	Special Instructions For Lab	g							
Notes: 1.) Fraction T = Tota 2.) QC Codes: FS = F 3.) Sample Markix G) 4.) Preservation Type 6.) Bottle Type: G = C	Notes: 1.) Fraction T = Total. D = Dissolved, S = SPLP, C = TCLP. N = Not Applicable 1.) Fraction T = Total. D = Dissolved, S = SPLP, C = TCLP. N = Not Applicable 2.) QC Codes: FS = Field Sample, TB = Trip Blank, FD = Field Duplicate, EB = Equipment Blank, MS = Matrix Spike, MSD = Matrix Spike, Duplicate, TB = Performance Evaluation Sample, FB = Field Blank 3.) Sample Matrix GW = Groundweter, SW = Surface Vester, DW = Driving Water SO = Soil, SD = Sedment, BW = Blank Water, NAL = Non-Aquacue Liquid, PR = Product O = Oil 4.) Preservation Type: Ha = Hydrochloric Acid, NI = Nitric Acid, SA = Sufficie Acid, SH = Soilum Hydroxide, Zn = Znr Acidate, ME = Merhanol, DI = DI Water 5.) Bottle Type: G = Glass P = Plastic, V = 4Gnt, VQA Glass Vial, AG = Amber Glass, AV = 4Gnt, VQA Amber Glass, Vial	C = TCLP. N = No. N, FD = Field Dup rface Visite. DW = = Nrtrc Acid SA , VOA Glass Visit.	N Applicable licate, E8 = Drinking V Sufuric A AG = Ambe	e Equipmen Nater SO ** Acid SH ** or Glass, A)	tt Blank M: Soil, SD ÷ Sodium 11y	S = Malrix E Sediment: droxide, Zn	ipka, MSD = BW = Blank = Znc Acela	Matrix Spiri Water, NAL	a Duylcala * Non-Aqu ethanol, Dl	, PE » Perfor secus Liquid, = Dj Water	mence Evak PR = Produc	uskion Semp ¤O∗Oil	. FB × Fiet	d Blank				
															Cooler 7	Cooler ? 60 / N MADEP Requirement	P Requiremen	드

Samples Iced 80 N	p receip	By: Date:	J128	9
	/ / Time:	/1/11 Time:	3 8250 z	
	Date: /	Date: 6	94: 867	
	Received:	Received: howy by Date: 6 /7 / 1) Time: By:	ighed Fed Ex P1 + mction#: 8673 8250 3810	
	_ Time: 15a	Time:F	ripped Fed 69	
		Date://	ON Caler Sh	
†	The has		•	
	Relinquished:	Relinquished	Ř7Z.	98¥Z

AINO70/1250628/6310728-31

	Г.	-,-				1 1		1		
	CPreservative Type (4)	Datie 1/De (2)			Comments (Special Instructions)	Chin I'm tor	Ant Volume		LUNCASTEC LAB	page 20f2
						· ;	-		-	
			(SW/S//3	יַ רכ	168 noM)		í	1	-	
≥	√		HMCU 'FIN) - -	Hyd-azina		1	1		
ಶ	9	, <u></u>			TEN NOM)		•	1		
ogram	<u> </u>	: 	ниао ,ни	M ,s	Hyd: azin	×	<u>i×</u>	<u> </u> ×		
Ö						-	1.	1	-	
ر ب ا	. 					-	<u>:</u>	+	•	
+	•						i	 	•	
<u></u>							Ì			
6	,		·			<u> </u>	· ·	<u> </u>		
\mathcal{N}	<u> -</u>	_			 -		<u>-</u> !	<u> </u>	-	
3					 -		,	1	•	
į						-	-	1	•	
:								-	•	
							<u> </u>	<u>.</u>		
	. >		(0089	a) are	Perchlora	,—	!	 		
÷	• >	(၁)	J9H - 80008) 910			;		!- -		
i i	· >		33 - GCMPD))8 bd	DMF (Mc	-		 		
1	· o		33 - CCW ₂ D)			=	: -	1		
	٠ م				17) 9+10		<u> </u>	<u> </u>		
į	· 0				Tot:#:40			<u>. </u>		
			C) or Grab (G)			-		<u> </u>		
ļ.			(3)	ieM i	San ple	<u> </u>	<u>:0</u>	FS SOI G		
4					उट द∞	Š	0 ∑ %	<u> </u>		
9			··			D CS W T	<u> ∑ ∨</u>	<u>18</u>		
_				1:11	ocapen4	<u>-</u>	<u> </u>	<u>i</u> -		
:)			Data:Time Collected	6/8/2011 9:35:00 AM	9 A	6/8/2011 9:35:00 AM		
į	I)			돌를	3/8/2 5.00	5.07	5.07		
		다. 교			20	9.69	<u>6</u>	.96		
;	ř -	י ק			<u> </u>	8	ည္			
7	MACTED	/ *1		i		OC-SS-448-0.0/1.0-MSD	OC-SS-448-0.071.0-XMS	OC-SS-448-0.0/1.0-XXX		
1	⋖	4				0,7	į	6		
	7	■ [_	9	(O)) (i) (ii)		
		•		}	0 4	S-44	244	4		
į					Sample ID	Ö	် ပ	S		
-			<u> </u>		<u></u>	<u>,Ó</u>	<u>.ŏ</u>	<u>10</u>		

Special Instructions For Lab	1) Fraction: T = Total D = Dissolved. S = SPLP, C = TCLP, N = Not Applicable	2.) OC Codes I FS = Field Sample, TB = Trip Blank, FD = Field Dupktorie (ED = Equipment Blank, MS = Maria Spike, MSD = Matta Spike Dupktorie, PE = Performance Evaluation Sample, FB = Field Blank	3.) Sample Matrix: GW = Groundwater, SW = Surface Water, DW = Drinking Water SO = Soil. SD = Sediment, BW = Benk Water, NAL = Non-Aqueous Liquid, PR = Product, O = Oil	4) Preservation Type: HA = Hydrochreric Acid, Ni = Nutric Acid, SA = Suffuric Acid, SH = Sociate Hydroxide, Zn = Zinc Acetate, ME = Nethanol, Di = Di Weter	
------------------------------	---	--	---	---	--

Cooler 20 / N MADEP Requirement	Samples loed (V) N	Temp @ receipt: 1-8 Deg C	Time Preservation / pH checked? Y / N	
		Time:	1 me	306
		Date:	Date: h / a / 1) Time: last	
		Received:	Received: H. / H.	
	!!	11 Time: 150 Received:	Time:	
	(Date: 6/8/	Date	
	F	Relinquished: Dis Mut	Relinquished:	E.A.

72 8583



Environmental Sample Administration Receipt Documentation Log

Client/l	Project: <u>0</u> '	lin Cup		Shippin	g Containe	er Sealed; (YĒ	S) NO				
Date of	Receipt: _	019111		Custody	y Seal Pres	sent*: YE	s) no				
Time o	f Receipt:	905		* Custody	seal was inta	ct unless otherwise					
Source	Code:	20-J		Packago	liscrepancy se B:	Chilled	Not Chilled				
			Temperature of	Shipping Contai	iners						
Cooler #	Thermometer ID	Temperature (*C)	Temp Bottle (TB) or Surface Temp (ST)	Wet Ice (WI) or Dry Ice (DI) or Ice Packs (IP)	Ice Present? Y/N	Loose (L) Bagged Ice (B) or NA	Comments				
1	ayes	1.8°L	TB	WI	ч	В					
2											
3					·						
4											
5	5										
6											
Number	r of Trlp Blank	s received <u>N</u>	OT listed on chain	of custody	0						
Paperw	ork Disc <i>r</i> epa	ncy/Unpack	ing Problems:				-				
					·						
		<u> </u>									
						·					
Unpack	er Signature <i>i</i>	/Emp#:	Many Hel	2316	Date/Tir	ne: <u>(alali\</u>	MS OFFICE BEBR				

Issued by Dept. 6042 Management

2174.06



Method Summary/Reference for SDG# OLN72_I

2425 New Holland Pike, PO Box 12425, Lancaster, PA 17605-2425 · 717-656-2300 Fax: 717-656-2681 · www.lancasterlabs.com

10346 Hydrazines in Soil

The soil is extracted with a buffer solution of known pH. An aliquot of the supernatant is derivatized and directly analyzed by HPLC/MS/MS.

Reference: Test Methods for Evaluating Solid Wastes, SW-846 Method 8315A modified, December 1996.

00111 Moisture

00118 Moisture

00121 Moisture Duplicate

A well-mixed sample is placed in a tared container and dried to a constant weight in an oven at 103-105C. The increase in weight is the total solids.

Reference: Standard Methods for the Examination of Water and Wastewater, $20 \, \text{th}$ Edition, 1998, Method $2540 \, \text{G}$

2425 New Holland Pike, PO Box 12425, Lancaster, PA 17605-2425 -717-656-2300 Fax: 717-656-2681 - www.kancasterlabs.com

ANALYTICAL RESULTS

Prepared by:

Prepared for:

Lancaster Laboratories 2425 New Holland Pike Lancaster, PA 17605-2425

Olin Corporation Suite 200 3855 North Ocoee Street Cleveland TN 37312

June 17, 2011

Project: Olin Wilmington, MA Superfund Site/6107090016

Submittal Date: 06/09/2011 Group Number: 1250628 SDG: OLN72 PO Number: REWI0012 Release Number: ERRE9813 State of Sample Origin: MA

Client Sample Description	Lançaster Labs (LLI) #
OC-SS-448-0.0/1.0-DUP Grab Soil	6310728
OC-SS-448-0.0/I.0-XXX Grab Soil	6310729
OC-SS-448-0.0/1.0-XMS Grab Soil	6310730
OC-SS-448-0.0/1.0-MSD Grab Soil	6310731

The specific methodologies used in obtaining the enclosed analytical results are indicated on the Laboratory Sample Analysis Record.

ELECTRONIC

MACTEC

Attn: Kelly Chatterton

COPY TO

ELECTRONIC

Attn: Chris Ricardi

COPY TO

MACTEC

ELECTRONIC

Olin Chemicals

Attn: James Cashwell

COPY TO

ELECTRONIC

Data Package Group

COPY TO



2425 New Holland Pike, PO Box 12425, Lancaster, PA 17605-2425 *717-656-2300 Fex: 717-656-2681 * www.lancasterlabs.com

Questions? Contact your Client Services Representative Nicole L Maljovec at (717) 656-2300 Ext. 1537

Respectfully Submitted,

Dorothy M. Love Group Leader

Douthy M. dove



Explanation of Symbols and Abbreviations

The following defines common symbols and abbreviations used in reporting technical data:

RL N.D. TNTC IU	Reporting Limit none detected Too Numerous To Count International Units	BMQL MPN CP Units NTU	Below Minimum Quantitation Level Most Probable Number cobalt-chloroplatinate units nephelometric turbidity units
umhos/cm	micromhos/cm	ng	nanogram(s)
· с	degrees Celsius	F	degrees Fahrenheit
meq	milliequivalents	lb.	pound(s)
g	gram(s)	kg	kilogram(s)
ug	microgram(s)	mg	milligram(s)
ml	milliliter(s)	Ī	liter(s)
m3	cubic meter(s)	ul	microliter(s)

- less than The number following the sign is the <u>limit of quantitation</u>, the smallest amount of analyte which can be reliably determined using this specific test.
- > greater than
- J estimated value The result is ≥ the Method Detection Limit (MDL) and < the Limit of Quantitation (LOQ).
- ppm parts per million One ppm is equivalent to one milligram per kilogram (mg/kg), or one gram per million grams. For aqueous liquids, ppm is usually taken to be equivalent to milligrams per liter (mg/l), because one liter of water has a weight very close to a kilogram. For gases or vapors, one ppm is equivalent to one microliter of gas per liter of gas.
- ppb parts per billion
- Dry weight Results basis concern

Results printed under this heading have been adjusted for moisture content. This increases the analyte weight concentration to approximate the value present in a similar sample without moisture. All other results are reported on an as-received basis.

Inorganic Qualifiers

U.S. EPA CLP Data Qualifiers:

Organic Qualifiers

Α	TIC is a possible aldol-condensation product _	В	Value is <crdl, but="" th="" ≥idl<=""></crdl,>
В	Analyte was also detected in the blank	E	Estimated due to interference
C	Pesticide result confirmed by GC/MS	M	Duplicate injection precision not met
D	Compound quantitated on a diluted sample	N	Spike sample not within control limits
. E	Concentration exceeds the calibration range of	S	Method of standard additions (MSA) used
	the instrument		for calculation
N	Presumptive evidence of a compound (TICs only)	U	Compound was not detected
· P	Concentration difference between primary and	W	Post digestion spike out of control limits
	confirmation columns >25%	*	Duplicate analysis not within control limits
U	Compound was not detected	+	Correlation coefficient for MSA < 0.995
Y,Z	Defined in case narrative		

Analytical test results meet all requirements of NELAC unless otherwise noted under the individual analysis.

Measurement uncertainty values, as applicable, are available upon request.

Tests results relate only to the sample tested. Clients should be aware that a critical step in a chemical or microbiological analysis is the collection of the sample. Unless the sample analyzed is truly representative of the bulk of material involved, the test results will be meaningless. If you have questions regarding the proper techniques of collecting samples, please contact us. We cannot be held responsible for sample integrity, however, unless sampling has been performed by a member of our staff. This report shall not be reproduced except in full, without the written approval of the laboratory.

WARRANTY AND LIMITS OF LIABILITY - In accepting analytical work, we warrant the accuracy of test results for the sample as submitted. THE FOREGOING EXPRESS WARRANTY IS EXCLUSIVE AND IS GIVEN IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED. WE DISCLAIM ANY OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING A WARRANTY OF FITNESS FOR PARTICULAR PURPOSE AND WARRANTY OF MERCHANTABILITY. IN NO EVENT SHALL LANCASTER LABORATORIES BE LIABLE FOR INDIRECT, SPECIAL, CONSEQUENTIAL, OR INCIDENTAL DAMAGES INCLUDING, BUT NOT LIMITED TO, DAMAGES FOR LOSS OF PROFIT OR GOODWILL REGARDLESS OF (A) THE NEGLIGENCE (EITHER SOLE OR CONCURRENT) OF LANCASTER LABORATORIES AND (B) WHETHER LANCASTER LABORATORIES HAS BEEN INFORMED OF THE POSSIBILITY OF SUCH DAMAGES. We accept no legal responsibility for the purposes for which the client uses the test results. No purchase order or other order for work shall be accepted by Lancaster Laboratories which includes any conditions that vary from the Standard Terms and Conditions, and Lancaster hereby objects to any conflicting terms contained in any acceptance or order submitted by client.



2425 New Holland Pike, PO Box 12425, Lancaster, PA 17605-2425 •717-656-2300 Fax: 717-656-2681 • www.lancasterlabs.com

Page 1 of 1

Sample Description: OC-SS-448-0.0/1.0-DUP Grab Soil

Wilmington MA Superfund Site

LLI Sample # SW 6310728 LLI Group # 1250628

Account # 12506

Project Name: Olin Wilmington, MA Superfund Site/6107090016

Collected: 06/08/2011 09:35

Submitted: 06/09/2011 09:05

SDG#: OLN72-01FD

Olin Corporation

Suite 200

3855 North Ocoee Street

Cleveland TN 37312

Reported: 06/17/2011 14:11

448-D

CAT No.	Analysis Name		CAS Number	Dry Result	Dry Limit of Quantitation*	Dry Method Detection Limit	Dilution Factor
Misc.	Organics	SW-846 831	5 A	ng/g	ng/g	ng/g	
	_	modified					
10346	1.1-Dimethylhydrazi	ne	57-14- 7	N.D.	5.3	. 2.1	1
10346	Hydrazine		302-01-2	2.5	2.1	0.53	1
10346	Methylhydrazine		60-34-4	N.D.	5.3	2.1	1
Wet C	hemistry	SM20 2540	G	*	%	*	
00111	Moisture		n.a.	6.2	0.50	0.50	1
	"Moisture" represen 103 - 105 degrees C as-received basis.						

General Sample Comments

All QC is compliant unless otherwise noted. Please refer to the Quality Control Summary for overall QC performance data and associated samples.

Laboratory Sample Analysis Record

CAT No.	Analysis Name	Method	Trial#	Batch#	Analysis Date and Tim	18	Analyst	Dilution Factor
	Hydrazines in Soil	SW-846 8315A modified	1	11165002	06/16/2011	00:03	Meng Yu	1
00111	Moisture	SM20 2540 G	2	11165820002B	06/14/2011	19:21	Scott W Freisher	1



2425 New Holland Pike, PO Box 12425, Lancaster, PA 17605-2425 •717-656-2300 Fax: 717-656-2681 • www.lancasterlabs.com

Page 1 of 1

Sample Description: OC-SS-448-0.0/1.0-XXX Grab Soil

Wilmington MA Superfund Site

LLI Sample # SW 6310729 LLI Group # 1250628

Account # 12670

Project Name: Olin Wilmington, MA Superfund Site/6107090016

Collected: 06/08/2011 09:35

Reported: 06/17/2011 14:11

Olin Corporation

Suite 200

Submitted: 06/09/2011 09:05

3855 North Ocoee Street

Cleveland TN 37312

448-1 SDG#: OLN72-02BKG

CAT No.	Analysis Name		CAS Number	Dry Result	Dry Limit of Quantitation*	Dry Method Detection Limit	Dilution Factor
Misc.	Organics	5W-846 83	15A	ng/g	ng/g	ng/g	
		modified					
10346	1,1-Dimethylhydra	zine	57-14-7	N.D.	5.3	2.1	1
10346	Hydrazine		302-01-2	1.7 J	2.1	0.53	1
10346	Methylhydrazine		60-34-4	N.D.	5.3	2.1	1
Wet C	hemistry	SM20 2540	G	*	*	%	
00111	Moisture .		n.a.	6.5	0.50	0.50	1
	"Moisture" repres 103 - 105 degrees as-received basis	Celsius. The r					

General Sample Comments

All QC is compliant unless otherwise noted. Please refer to the Quality Control Summary for overall QC performance data and associated samples.

Laboratory Sample Analysis Record

CAT No.	Analysis Name	Method	Trial#	Batch#	Analysis Date and Time	Analyst	Dilution Factor
10346	Hydrazines in Soil	SW-846 8315A modified	1	11165002	06/15/2011 23:43	Meng Yu	1
00111	Moisture	SM20 2540 G	2	11165820002B	06/14/2011 19:21	Scott W Freisher	1



2425 New Hotland Pike, PO Box 12425, Lancaster, PA 17605-2425 • 717-656-2300 Fax: 717-656-2681 • www.lancasterlabs.com

Page 1 of I

Sample Description: OC-SS-448-0.0/1.0-XMS Grab Soil

Wilmington MA Superfund Site

LLI Sample # SW 6310730 LLI Group # 1250628

Account # 12670

Project Name: Olin Wilmington, MA Superfund Site/6107090016

Collected: 06/08/2011 09:35

Olin Corporation

Suite 200

Submitted: 06/09/2011 09:05 3855 North Ocoee Street

Reported: 06/17/2011 14:11 Cleveland TN 37312

448-1 SDG#: OLN72-02MS

CAT No.	Analysis Name	CAS Number	Dry Result	Dry Limit of Quantitation*	bry Method Detection Limit	Dilution Factor
Misc.	Organics	SW-846 8315A modified	ng/g	ng/g	ng/g	
10346	1,1-Dimethylhydrazi	ne 57-14-7	59	5.3	2.1	1
10346	Hydrazine	302-01-2	15	2.1	0.53	1
10346	Methylhydrazine	60-34-4	54	5.3	2.1	1
Wet C	hemistry	SM20 2540 G	%	*	*	
00118	Moisture	n.a	6.5	0.50	0.50	1

General Sample Comments

All QC is compliant unless otherwise noted. Please refer to the Quality Control Summary for overall QC performance data and associated samples.

Laboratory	Samnle	Analveie	Pecord
Laboratory	Samble	WIIGTARTR	RECOLG

CAT No.	Analysis Name	Method	Trial#	Batch#	Analysis Date and Tim	ie	Analyst	Dilution Factor
10346	Hydrazines in Soil	SW-846 B315A modified	1	11165002	06/16/2011	01:45	Meng Yu	1
00118	Moisture	SM20 2540 G	2	11165820002B	06/14/2011	19:21	Scott W Freisher	1



2425 New Holland Pike, PO Box 12425, Lancaster, PA 17605-2425 •717-556-2300 Fax: 717-556-2681 • www.lancasterlabs.com

Page 1 of 1

Sample Description: OC-SS-448-0.0/1.0-MSD Grab Soil

Wilmington MA Superfund Site

LLI Sample # SW 6310731 LLI Group # 1250628

Account # 12670

Project Name: Olin Wilmington, MA Superfund Site/6107090016

Collected: 06/08/2011 09:35

Olin Corporation

Suite 200

Submitted: 06/09/2011 09:05 3855 North Ocoee Street

Cleveland TN 37312

Reported: 06/17/2011 14:11

448-1 SDG#: OLN72-02MSD*

CAT No.	Analysis Name		CAS Number	Dry Result	Dry Limit of Quantitation*	Dry Method Detection Limit	Dilution Factor
Misc.	Organics	sw-846 831	.5 A	ng/g	ng/g	ng/g	
	_	modified					
10346	1,1-Dimethylhydrazin	e	57-14-7	60	5.3	2.1	1
10346	Hydrazine		302-01-2	13	2.1	0.53	1
10346	Methylhydrazine		60-34-4	39	5.3	2.1	1
Wet C	hemistry	SM20 2540	G	*	4	•	
00118	Moisture		n.a.	6.5	0.50	0.50	1
00121	Moisture Duplicate		n.a.	6.1	0.50	0.50	1
	The duplicate moistumoisture test. For determination is the	comparability	y purposes, ti	he initial mod	isture		

General Sample Comments

All QC is compliant unless otherwise noted. Please refer to the Quality Control Summary for overall QC performance data and associated samples.

Laboratory Sample Analysis Record

CAT No.	Analysis Name	Method	Trial#	Batch#	Analysis Date and Tim	10	Analyst	Dilution Factor
10346	Hydrazines in Soil	SW-846 8315A modified	1	11165002	06/16/2011	02:05	Meng Yu	1
00118	Moisture	SM20 2540 G	2	11165820002B	06/14/2011	19:21	Scott W Freisher	1
00121	Moisture Duplicate	SM20 2540 G	2	11165820002B	06/14/2011	19:21	Scott W Freisher	1

Hydrazines by LC/MS/MS

0

Case Narrative Conformance/Nonconformance Summary



CLIENT: Olin Corporation

SDG: OLN72

Specialty Services Group Fraction: Hydrazines by LC/MS/MS

Hydrazines in Soil

Ma	trix

Sample #	Client ID	<u>Liquid</u>	Solid	<u>Comments</u>
6310728	OC-SS-448-0.0/1.0-DUP		X	Field Duplicate Sample
6310729	OC-SS-448-0.0/1.0-XXX		X	Unspiked
6310730	OC-SS-448-0.0/1.0-XMS		X	Matrix Spike
6310731	OC-SS-448-0.0/1.0-MSD		X	Matrix Spike Duplicate

See QC Reference List for Associated Batch QC Samples

SAMPLE PREPARATION:

Samples were derivatized with benzaldehyde prior to analysis.

ANALYSIS:

There were no dilutions performed for analyses associated with samples in this SDG.

No problems were encountered with the analysis of the samples.

QUALITY CONTROL AND NONCONFORMANCE SUMMARY:

Please note that US EPA Methods for organic compounds do not require action by the laboratory based on out-of-specification MS/MSD results.

DATA INTERPRETATION:

No further interpretation is necessary for the data submitted.

DEATZ BELS

7/1/2011 9:59:35 AM Page 1 of 2



CLIENT: Olin Corporation SDG: OLN72

Specialty Services Group Fraction: Hydrazines by LC/MS/MS

Abbreviation Key

110,01011111111111111111111111111111111	
UNSPK = Unspiked (for MS/MSD)	LOQ = Limit of Quantitation
MS = Matrix Spike	MDL = Method Detection Limit
MSD = Matrix Spike Duplicate	ND = Not Detected
BKG = Background (for Duplicate)	J = Estimated Value
D = Duplicate (DUP)	E= out of calibration range
LCS = Lab Control Sample	
LCSD = Lab Control Sample Duplicate	* = Out of Specification
NC = Not calculated	NF = Not found

Narrative Reviewed and Approved 7/1/11 by Dowth M dwe.
(Date)

QC Summary



Quality Control Reference List Specialty Services Group

CLIENT: Olin Corporation

SDG: OLN72

Fraction: Hydrazines by LC/MS/MS

Analysis		
Hydrazines	in	Soil

Batch	Number
1116500	02

Sample Number	Analysis Date
BLK	06/15/2011 22:22:00
LCS	06/16/2011 01:04:00
LCSD	06/16/2011 01:25:00
6310728	06/16/2011 00:03:00
6310729 UNSPK	06/15/2011 23:43:00
6310730 MS	06/16/2011 01:45:00
6310731 MSD	06/16/2011 02:05:00



Fraction: Hydrazines by LC/MS/MS

Quality Control Summary Method Blank Specialty Services Group SDG: OLN72 Matrix: SOLID

11165002 / BLK Analyte	Analysis Date	Blank Results	Units	MDL	LOQ
Hydrazine	06/15/11	N.D.	ng/g	0.50	2.0
Methylhydrazine	06/15/11	N.D.	ng/g	2.0	5.0
1,1-Dimethylhydrazine	06/15/11	N.D.	ng/g	2.0	5.0



Quality Control Summary Matrix Spike/Matrix Spike Duplicate

SDG: OLN72 Matrix: SOLID

Specialty Services Group Fraction: Hydrazines by LC/MS/MS

UNSPK: 6310729	Batch: 111	65002 (Sampl	e number(s):	6310728-63	310731)				
MS: 6310730 MSD: 6310731 Analyte	Spike Added ng/g	Unspiked Conc ng/g	MS Conc ng/g	MSD Conc ng/g	MS %Rec	MSD %Rec	%Rec Limits	%RPD	%RPD Limits
Hydrazine	120	1.58	13.77	12.56	10 *	9 *	11-102	9	30
Methylhydrazine	600	N.D.	50.61	36.38	8 *	6 *	10-92	33 *	30
1,1-Dimethylhydrazine	600	N.D.	55.45	55.84	9 *	9 * _	10-116	1	30



Quality Control Summary Laboratory Control Standard (LCS) Laboratory Control Standard Duplicate(LCSD)

SDG: OLN72 Matrix: SOLID

Specialty Services Group Fraction: Hydrazines by LC/MS/MS

LCS	Batch: 111650	02 (Sample nur	nber(s): 631072	8-6310731)			
LCSD	Spike	LCS	LCSD					
	Added	Conc	Conc	LCS	LCSD	%Rec		%RPD
Analyte	ng/g	ng/g	ng/g	%Rec	%Rec	Limits	%RPD	Limits
Hydrazine	120	92.06	113.59	77	95	70-130	21	30
Methylhydrazine	600	424.3	481.02	71	80	70-130	13	30
1,1-Dimethylhydrazine	600	516.92	548.7	86	91	70-130	6	30

Sample Data



Fraction: Hydrazines by LC/MS/MS

LOQ/MDL Summary Specialty Services Group

SDG: OLN72

10346: Hydrazines in Soil Analyte Name	Default MDL	Default LOQ	Units_
Hydrazine	0.50	2.0	ng/g
Methylhydrazine	2.0	5.0	ng/g
1,1-Dimethylhydrazine	2.0	5.0	ng/g



Component Name:

Monomethylhydrazine

i	!	S	Summary of Quan Results	an Results				ĺ
Sample ID	Data File Name	Arca	ISTD Area	Area Ratio Specified Amount	mount	Calculated Amount	% Diff	Excluded
conditioner	A11165002 01	N/A	N/A	N/A	N/A	N/A	N/A	V/N
conditioner	A11165002 02	N/A	N/A	N/A	N/A	N/A	A/N	Υ/Z
SYS(MDL)	A11165002_03	24970.60	N/A	24970.595	N/A	2.77233ug/kg	A/N	V/Z
CAL1	A11165002_04	58896,46	A/N	58896.458	S	6.40898ug/kg	28.18	N/A
CAL2	A11165002 05	105001.17	N/A	105001.171	10	11.35113ug/kg	13.51	V/N
CAL3	A11165002_06	183559.44	V/N	183559.444	22	19.77211ug/kg	-20.91	A/A
CAL4	A11165002_07	415437.59	N/A	415437.590	50	44.62807ug/kg	-10.74	Y/V
CAL5	A11165002_08	2072605.08	N/A	2072605.077	250	222.26658ug/kg	-11.09	Y/A
CAL6	A11165002 09	4504468.74	N/A	4504468.743	200	482.94791ug/kg	-3.41	N/A
CAL7	A11165002_10	9448410.12	N/A	9448410.120	1000	1012.90905ug/kg	1.29	N/A
CAL8	A11165002_11	12030709.14	N/A	12030709.138	1250	1289.71616ug/kg	3.18	A/N
Conditioner	A11165002 12	N/A	N/A	Ϋ́Х	N/A	Y/Z	A/N	A/N
Conditioner	A11165002 13	N/A	V/N	N/A	N/A	V/Z	V/V	N/A
BLK Sand	A11165002 14	N/A	Y/X	N/A	0	N/A	Y/A	Y/A
LACO	A11165002_15	204017.63	N/A	204017.632	53	21.96511ug/kg	-12.14	V/V
Conditioner	A11165002 16	N/A	N/A	N/A	N/A	V/A	N/A	V/A
6310729(BKG)	A11165002 ⁻ 17	N/A	N/A	N/A	N/A	N/A	N/A	Y/V
6310728	A11165002_18	N/A	N/A	N/A	N/A	N/A	N/A	V/V
CCV2	A11165002 19	398473.08	V/N	398473.082	8	42.80958ug/kg	-14.38	N/A
ICV	A11165002_20	4875265.14	A/N	4875265.142	∀ Z	522.69508ug/kg	N/A	Y/X
TCS	A11165002_21	3957319.27	Y/N	3957319.268	Ν	424.29674ug/kg	N/A	Y/Z
TCSD	A11165002_22	4486479.01	N/A	4486479.012	Y/Z	481.01952ug/kg	V/V	N/A
6310730 (MS)	A11165002 23	471272.24	N/A	471272.236	N/A	50.61322ug/kg	N/A	V/V
6310731(MSD)	A11165002 24	338527.47	N/A	338527.471	N/A	36.38377ug/kg	N/A	N/A
CCV3	A11165002_25	2076710.46	N/A	2076710.458	250	222.70665ug/kg	-10.92	V/Z
V/N	A11165002_26	4588763.13	N/A	4588763.128	200	491.98377ug/kg	-1.60	Y /Z

Page 1 of 3 Thursday, June 16, 2011, 12:53:54



Component Name:

1,1-Dimethylhydrazine

		5 21	Summary of Quan Results	an Results				
Sample ID	Data File Name	Arca	ISTD Area	Area Ratio Sper	Specified Amount	Catculated Amount	% Diff	Excluded
conditioner	A11165002 01	N/A	N/A	N/A	N/A	N/A	N/A	V/N
conditioner	A11165002 02	N/A	N/A	N/A	N/A	N/A	N/A	Y/X
SYS(MDL)	A11165002 03	24427.88	N/A	24427.879	N/A	3.83821ug/kg	N/A	N/A
CALI	A11165002 04	50181.01	N/A	50181.009	S	6.44055ug/kg	28.81	N/A
CAL2	A11165002_05	112117.31	N/A	112117,312	91	12.69919ug/kg	26.99	N/A
CAL3	A11165002 06	167545.26	N/A	167545.255	25	18.30016ug/kg	-26.80	N/A
CAL4	A11165002_07	381544.38	Y/X	381544.385	50	39,92469ug/kg	-20.15	N/A
CALS	A11165002_08	2112373.87	Y/X	2112373,869	250	214.82433ug/kg	-14.07	N/A
CAL6	A11165002 09	5003884.12	N/A	5003884.123	200	507.01031ug/kg	1.40	Y/N
CAL7	A11165002_10	10155323.34	N/A	10155323.337	1000	1027.56122ug/kg	2.76	A/N
CAL8	A11165002_11	12487626.54	N/A	12487626.542	1250	1263.23955ug/kg	90'1	A/N
Conditioner	A11165002_12	11189.34	V/Z	11189.337	N/A	2.50046ug/kg	N/A	N/A
Conditioner	A11165002_13	N/A	A/N	A/N	A/N	Y/X	N/A	N/A
BLK Sand	A11165002_14	N/A	A/N	N/A	0	V/N	N/A	N/A
CCVI	A11165002 15	.169966.34	V/N	169966.337	25	18.54481ug/kg	-25.82	A/N
Conditioner	A11165002 16	N/A	A/N	N/A	N/A	N/A	A/A	N/A
6310729(BKG)	A11165002_17	N/A	A/N	N/A	N/A	Y / X	N/A	Y/A
6310728	A11165002_18	N/A	N/A	N/A	N/A	N/A	A/N	N/A
CCV2	A11165002_19	402245.48	V/N	402245.481	20	42.01653ug/kg	-15.97	N/A
ICV	A11165002 20	5614302.18	N/A	5614302.178	N/A	568.69281ug/kg	K/Z	Y/X
rcs	A11165002_21	5101969.82	N/A	\$101969.819	N/A	516.92183ug/kg	K/N	A/A
CSD	A11165002_22	5416398.55	A/N	5416398,549	N/A	548.69473ug/kg	A/N	A/Z
6310730 (MS)	A11165002 23	535146.49	A/N	535146.490	A/N	55.44612ug/kg	N/A	V/Z
6310731(MSD)	A11165002 24	539001.85	A/N	539001.853	N/A	55.83570ug/kg	N/A	ΥX
CCV3	A11165002_25	2163973.65	ΥX	2163973,647	250	220.03847ug/kg	-11.98	A/Z
N/A	A11165002_26	4940768.03	N/A	4940768.035	200	500.63245ug/kg	0,13	N/A

Page 2 of 3 Thursday, June 16, 2011, 11:36:00



Hydrazine

Component Name:

LCMSMS ANALYSIS REPORT

% Diff Excluded				16.98 N/A	-3.36 N/A		-1.81 N/A						N/A N/A		5.99 N/A							N/A N/A		N/A N/A		
Calculated Amount	N/A	N/A	0.74426ug/kg	1.16975ug/kg	1,93289ug/kg	4.55090ug/kg	9.81926ug/kg	48.77498ug/kg	98.35341ug/kg	199.06698ug/kg	254.33183ug/kg	A/N	V /N	V/V	5.29948ug/kg	N/A	1.57732ug/kg	2.36091ug/kg	9.26236ug/kg	121.91427ug/kg	92.05647ug/kg	113.59169ug/kg	13.76628ug/kg	12.56267ug/kg	50.42795ug/kg	101.95675ug/kg
Specified Amount	N/A	N/A	N/A	_	7	Š	10	20	100	200	250	N/A	N/A	0	\$	N/A	N/A	N/A	10	N/A	N/A	N/A	N/A	N/A	50	100
tio	N/A	V/N	8361.611	16210.294	30287,338	78579.410	175760.585	894343.943	1808875.647	366654.195	4686078.544	N/A	V/N	V/V	92387.898	N/A	23728.397	.38182.657	165487.879	2243482.946	1692721.302	2089963.356	248567.911	226365,897	924834.870	1875343.421
ISTD Area Area Ra	N/A	N/A	A/N	A/X	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	A/N	Ϋ́Χ	V/N	A/N	N/A	V /N	N/A	N/A	N/A	A'N	V /N	N/A	A/N	Y/N
Area	N/A	N/A	8361.61	16210.29	30287.34	78579.41	175760.58	894343.94	1808875.65	3666654.19	4686078.54	N/A	A/A	Y.Z	92387.90	A/N	23728.40	38182.66	165487.88	2243482.95	1692721.30	2089963.36	248567.91	226365.90	924834.87	1875343 42
Data File Name	A11165002 01	A11165002_02	A11165002_03	A11165002 04	A11165002 05	A11165002_06	A11165002 07	A11165002 08	A11165002 09	A11165002 10	A11165002 11	A11165002 12	A11165002 13	A11165002_14	A11165002_15	A11165002 16	A11165002 17	A11165002 18	A11165002_19	A11165002_20	A11165002_21	A11165002_22	A11165002 23	A11165002_24	A11165002_25	A11165002 26
Sample ID	conditioner	conditioner	SYS(MDL)	CALI	CAL	CAL3	CAL4	CALS	CAL6	CAL7	CAL8	Conditioner	Conditioner	BLK Sand	CCV1	Conditioner	6310729(BKG)	6310728	CCV2	ICV	CCS	CSD	6310730 (MS)	6310731(MSD)	CCV3	N/A

Page 3 of 3 Thursday, June 16, 2011, 11:36:00



Sample Name:

6310728

Data File:

A11165002_18

Sample Type: Run Time(min): Unknown

Injection Volume(µI):

9.99

Dilution Factor: Instrument Model:

Instrument Method:

5.00 1.00

TSQ Quantum Access C:\XCalibur\Hydrazine

Analysis\Hydraz_soil

Acquisition Date:

Sample ID:

Vial:

06/16/11 12:03:58 AM 6310728

a:19

Instrument Software Version: 1.4.1

Instrument Name:

Original Data Path:

Instrument Serial Number:

Ouantum

TOU01408

C:\XCalibur\Hydrazine


Analysis\2011June

Operator:


Quantum

Quan Peak Table


	V MARK W	THE PARTY NAMED IN COLUMN TO THE PARTY NAMED		
Component Name	Calculated Amount	Units	Response Ratio	 RT
1,1-Dimethylhydrazine	N/A	ug/kg	N/A	N/A
Monomethylhydrazine	N/A	ug/kg	N/A	N/A
Hydrazine	2.361	ug/kg	38182.657	8.46




There's no data available to display this graphic object.



There's no data available to display this graphic object.





Page 1 of 1 Thursday, June 16, 2011, 11:37:22



Sample Name:

6310729(BKG)

Data File:

A11165002 17

Sample Type:

Unknown

Run Time(min):

9.98

Injection Volume(µI):

5.00

Dilution Factor:

Instrument Model: Instrument Method: 1.00

TSO Quantum Access C:\XCalibur\Hydrazine

Analysis\Hydraz soil

Acquisition Date:

Sample ID:

06/15/11 11:43:41 PM

6310729(BKG)

a:18 1.4.1

Instrument Software Version:

Instrument Name:

Quantum

Instrument Serial Number: Original Data Path:

TQU01408

C:\XCalibur\Hydrazine

Analysis\2011June


Operator:

Quantum


Quan Peak Table

Vial:


Component Name	Calculated Amount	Units	Response Ratio	 RT
1,1-Dimethylhydrazine	N/A	ug/kg	N/A	N/A
Monomethylhydrazine	N/A	ug/kg	N/A	N/A
Hydrazine	1.577	ug/kg	23728.397	8.43



There's no data available to display this graphic object.



There's no data available to display this graphic object.



A11165002_17 #617_RT; 8.43_AV; 1_NL; 3.18E3 F; + c APCI SRM ms2 209.070@cid20.00 [77.325-77.335, 106.215-106.225] 77,33 100.220万工程7.206.28日2日 77.335 77.330

Page I of 1 Thursday, June 16, 2011, 11:37:21

Standards Data



						S ₂	Sequence Table		
File Name	Sample ID	Sample Type	Level	Vial	II S	Dil Factor	Path	Inst Method	Proc Method
A11165002_01	conditioner	Unknown	N/A	A:1	5.0 	1.000	CAXCaliburHydrazine Analysis/2011June	C:\XCalibur\Hydrazine Analysis\Hydraz_soil	· C:\XCalibur\Hydrazine Analysis\Processing Methods\Hydraz soil
A11165002_02	conditioner	Unknown	N/A	A:1	5.0	1.000	C:\XCalibur\Hydrazine Analysis\2011June	C:\XCalibur\Hydrazine Analysis\Hydraz_soil	C:\XCalbur\Hydraz_son C:\XCalbur\Hydrazinc Analysis\Processing Achode\Hydraz evil
A11165002_03	SYS(MDL)	Unknown	N/A	A:2	5.0	1.000	C:\XCalibur\Hydrazine Analysis\2011June	C:XCaliburMydarazine AnalysisMydraz_soil	C:\XCalibu\I;\ydrazinc C:\XCalibu\I;\ydrazinc Analys\Processing Methode\Fydraz soil
A11165002_04	CALI	Std Bracket	,- 1	A:3	5.0	1.000	C:\XCalibur\Hydrazine Analysis\2011June	C:\XCalibur\Hydrazinc Analysis\Hydraz_soil	C:\Caliburklydrazine C:\Analysis\Processing Methods\Pydraz.soil
A11165002_05	CAL2	Std Bracket	7	A:4	5.0	1.000	C:\XCalibur\Hydrazine Analysis\2011une	C:\XCalibur\Hydrazinc Analysis\Ffydraz_soil	C:\CCalibur\!gdrazinc Analysis\Processing Methods\Hydraz soil
A11165002_06	CAL3	Std Bracket	e	A:5	5.0	1.000	CAXCaliburHydrazine Analysis/2011June	C:\XCalibur\Hydrazine Analysis\Hydraz_soil .	C:\XCalibur\Iydrazine Analysis\Processing Methods\Hydraz soil
A11165002_07	CAL4	Std Bracket	4	A :6	5.0	1,000	C:\XCalibur\Hydrazinc Analysis\2011June	C:\XCalibur\Hydrazine Analysis\Hydraz_soil	C:\XCaliburkIydrazine Analysis\Processing Methods\Pydraz soil
A11165002_08	CALS	Std Bracket	5	A:7	5.0	1.000	C:\XCaliburHydrazine Analysis\2011June	C:\XCalibur\Hydrazine Analysis\Hydraz_soil	C:\XCalibur\Hydrazine Analysis\Processing Methods\Nutraz soil
A11165002_09	CAI.6	Std Bracket	9	A:8	5.0	1.000	CAXCalibur/Hydrazine Analysis/2011June	C:\XCalibur\Iydrazine Analysis\Iydraz_soil	C:\XCalibury\xxiazinc Analysis\Processing Methods\lvatez
A11165002_10	CAL7	Std Bracket	7	A:9	5.0	1,000	CAXCalibuvHydrazine Analysis/2011June	C:\XCalibur\Hydrazine Analysis\Hydraz_soil	C:\ZCalibux !ydrazine Analysis\Processing Methods\Hydraz soil
A11165002_11	CAL8	Std Bracket	∞	A:10	5.0	000'1	C:\XCaliburHydrazine Analysis\2011June	C:\XCalibur\Hydrazine Analysis\Hydraz_soil	C:\XCaliburHydrazine Analysis\Processing Achods\Ivdaz soil
A11165002_12	Conditioner	Unknown	N/A	a:1	5.0	1,000	CAXCaliburHydrazine Analysis/2011June	C:\XCalibur\Hydrazine Analysis\Hydraz_soil	C:\XCaliburHydrazine Analysis\Processing Archods\Hydraz soil
A11165002_13	Conditioner	Unknown	N/A	a :1	5,0	1.000	C:XCaliburHydrazine Analysis/2011June	C:\XCalibur\Hydrazine Analysis\Hydraz_soil	C:\Calcalogous and Control of Con
A11165002_14	BLK Sand	Blank	N/A	a:11	5.0	1.000	C:\XCalibur\Hydrazine Analysis\2011June	C:\XCalibur\Hydrazine Analysis\Hydraz_soil	C:XCallbary
8838								Page 1 of 2 Friday, July 01, 2011, 14:45:51	

		2							
File Name	Sample ID	Sample Type	Level	Vial	lnj Vol	Dil Factor	Path	Inst Method	Proc Method
A11165002_15	CCVI	90	_	a:5	5.0	1,000	C:\XCalibur\Hydrazine Analysis\2011June	C:\XCalibur\Hydrazine Analysis\Hydraz_soil	C:\XCalibut\Hydrazine Analysis\Processing Mothod@Hydraz soil
A11165002_16	Conditioner	Unknown	N/A	a: 1	5.0	1,000	C:\XCalibur\Hydrazine Analysis\2011June	C:\XCalibur\Hydrazine Analysis\Hydraz_soil	C./XCalibur/Hydrazine Analysis/Processing
A11165002_17	6310729(BKG)	Unknown	N/A	a:18	5.0	1.000	C:\XCalibur\Hydrazine Analysis\2011June	C:/XCalibu/Hydrazinc Analysis/Hydraz_soil	C:XCaliburIyanaz sen C:XCaliburIyadazine AnalysisVrocessing Matecdatt dense sedi
A11165002_18	6310728	Unknown	K/X	a:19	5.0	1.000	C:\XCalibur\Hydrazine Analysis\2011June	C:\XCalibur\Hydrazinc Analysis\Hydraz_soil	Wethoushtydrid, son CAXCaliburHydrazine Analysis/Processing Machadelladara soil
A11165002_19	CCV2	δ,	2	a:6	5.0	1.000	CAXCalibur/Hydrazine Analysis/2011June	CAXCaliburHydrazine AnalysisHydraz_soil	CAXCaliburHydrazine Analysis/Processing Machadell Lafax
A11165002_20	ICV	OC OC	ICV	a:13	5.0	1.000	C:\XCalibur\Hydrazine Analysis\2011June	C:\XCaliburHydrazine Analysis\Hydraz_soil	C:\XCulibur\Hydrazine Analysis\Processing Mathematical Analysis\Processing Mathematical Analysis\Processing Analysis\Processin
A11165002_21	гсѕ	Unknown	N/A	a:14	5.0	1.000	C:\XCalibur\Iydrazine Analysis\2011June	C:\XCalibur\Hydrazine Analysis\Hydraz_soil	Methodskitydnaz_soli C:XCalibur/Hydrazine Analysis/Processing Analysis/Processing
A11165002_22	TCSD	Unknown	N/A	a:15	5.0	1.000	C:\XCalibur\Iydrazine Analysis\2011June	C:\XCalibur\Hydrazine Analysis\Hydraz_soil	C:XCaliburHydrazine Analysis/Processing
A11165002_23	6310730 (MS)	Unknown	N/A	a:16	5.0	1.000	C:\XCalibur\Hydrazine Analysis\2011June	C:\XCalibur\Hydrazinc Analysis\Hydraz_soil	CAXCaliburAlydrazine Analysis/Processing
A11165002_24	6310731(MSD)	Unknown	N/A	a:17	5.0	1.000	C:\XCalibur\Hydrazine Analysis\2011June	C:\XCalibur\Hydrazine Analysis\Hydraz_soil	CAXCaliburHydrazine Analysis/Processing Mathodallurian
A11165002_25	CCV3		ш	a:7	5.0	1.000	C:\XCalibur\Hydrazine Analysis\2011June	C:\XCalibur\Hydrazine Analysis\Hydraz_soil	Memodavi yaraz_son C:XCaliburki yarazine Analysis/Processing Analysis/Processing
A11165002_26	CCV4	90	4	A:8	5.0	1.000	C:\XCalibur\Hydrazine Analysis\2011June	C:\XCalibu\Hydrazinc Analysis\Hydraz_soil	CARCAIBUNAINGEZING Analysis/Processing Methods/Hydraz_soil


Page 2 of 2 Friday, July 01, 2011, 14:45:51



Component Name:

Monomethylhydrazine

Monomethylhydrazine Y = -4166.25+9932.87*X R^2 = 0.9979 W: 1/X



Identification Filter:	+ c APCI SRM ms2 135.15@cid20.00	Component Name: 1st Trace Type:	Monomethylhydrazine Base Peak
2nd Trace Type: Mass Range 2 (m/z):	[77.33-77.33, 104.14-104.15] N/A	Mass Range 1 (m/z): Wavelength Range 2 (nm):	N/A
Base Peak(BP): Retention Time Window (sec): RT Reference: Adjust Using:	30.00000 No N/A	Expected RT (min): View Width (min): Adjust Expected RT:	3.74000 2.50000 No
Detection Options		Peak Detection Algorithm: ICIS Peak Integration	ICIS
ICIS Smoothing Points: Area Noise Factor: ICIS Constrain Peak Width; ICIS Tailing Factor:	3 5 No N/A	Baseline Window: Peak Noise Factor: ICIS Peak Height (%):	75 10 N/A
ICIS Peak Detection ICIS Minimum Peak Height (S/N): ICIS Window %:	50.0 N/A	ICIS Identify By: ICIS Ion Ratio Confirmation: ICIS Qualifier Ion Coelution (min): ICIS Spectrum Thresholds	Nearest RT N/A N/A
ICIS Forward: ICIS Match:	N/A N/A	ICIS Reverse:	N/A
ICIS Advanced Parameters Minimum Peak Width: Area Tail Extension:	3 5	Noise Method: Multiplet Resolution: Area Scan Window: Calibration	Incos 10 0
Component Type:	Target Compound	%RSD Calculation Method: Internal Standard	Use calculated amounts
ISTD Amount:	N/A	ISTD Units: Target Compounds	N/A
ISTD: Origin: Calibration Curve: Number of Cal. Levels:	IgnoreOrigin Linear 8	Weighting: Response: Target Units: Number of QC Levels: Peak Purity Options	OneOverX Area ug/kg 5
Scan Threshold (mAU): Limit ScanRange (nm):	N/A N/A	Peak Coverage (%):	N/A 7/5

Page 1 of 6 Friday, July 01, 2011, 15:15:46



Component Cal Level 7	<u> </u>
Cal Level	Amount
1	5.000
2	10.000
3	25.000
4	50.000
5	250.000
6	500.000
7	1000.000
8	1250.000

 Component OC Level Table

 QC Level
 Amount

 1
 25.000

 2
 50.000

 3
 250.000

 4
 500.000

 ICV
 600.000

	ICV & CCV R	esult Table		
Sample ID	Data File Name	Calculated Amount	Area	% Diff
CALI	A11165002 04	6.42269ug/kg	59629.55	28.45
CAL2	A11165002 05	11.13659ug/kg	106452.04	11.37
CAL3	A11165002 06	19.50917ug/kg	189615.91	-21.96
CAL4	A11165002 07	43.06328ug/kg	423575.90	-13.87
CAL5	A11165002 08	230.61361ug/kg	2286489.53	-7.75
CAL6	A11165002 09	505.03852ug/kg	5012317.47	1.01
CAL7	A11165002 10	1041.32862ug/kg	10339219.10	4.13
CAL8	A11165002 11	1232.88752ug/kg	12241949.44	-1.37
CCV1	A11165002 15	21.76000ug/kg	211973.10	-12.96
CCV2	A11165002 19	42.58962ug/kg	418871.00	-14.82
CCV3	A11165002 25	229.99733ug/kg	2280368.08	-8.00
CCV4	A11165002 26	500.76816ug/kg	4969900.47	0.15
ICV	A11165002_20	541.21537ug/kg	5371657.56	-9.80

OLNYZ JEHO


Page 2 of 6 Friday, July 01, 2011, 15:15:46



Component Name:

1,1-Dimethylhydrazine

1,1-Dimethylhydrazine Y = -13396.8+9885.75*X R^2 = 0.9963 W: 1/X



Identification Filter:	+ c APCI SRM ms2 149.10@cid30.00	Component Name: 1st Trace Type:	1,1-Dimethylhydrazine Base Peak
riker:	[77.33-77.33, 106.22-106.22]	ist itace type.	
2nd Trace Type: Mass Range 2 (m/z): Base Peak(BP):	N/A 106	Mass Range 1 (m/z): Wavelength Range 2 (nm):	N/A
Retention Time Window (sec): RT Reference: Adjust Using:	30.00000 No N/A	Expected RT (min): View Width (min): Adjust Expected RT:	5.91000 2.50000 No
Detection Options		Peak Detection Algorithm: ICIS Peak Integration	ICIS
ICIS Smoothing Points: Area Noise Factor: ICIS Constrain Peak Width: ICIS Tailing Factor:	3 5 No N/A	Baseline Window: Peak Noise Factor: ICIS Peak Height (%):	75 10 N/A
ICIS Peak Detection ICIS Minimum Peak Height (S/N): ICIS Window %:	100.0 N/A	ICIS Identify By: ICIS Ion Ratio Confirmation: ICIS Qualifier Ion Coelution (min): ICIS Spectrum Thresholds	Nearest RT N/A N/A
ICIS Forward: ICIS Match:	N/A N/A	ICIS Reverse:	N/A
ICIS Advanced Parameters Minimum Peak Width: Area Tail Extension:	3 5	Noise Method: Multiplet Resolution: Area Scan Window: Calibration	incos 10 0
Component Type:	Target Compound	%RSD Calculation Method: Internal Standard	Use calculated amounts
ISTD Amount:	N/A	ISTD Units: Target Compounds	N/A
ISTD: Origin: Calibration Curve: Number of Cal. Levels:	IgnoreOrigin Linear 8	Weighting: Response: Target Units: Number of QC Levels: Peak Purity Options Peak Coverage (%):	OneOverX Area ug/kg 5 02 N/A
Scan Threshold (mAU): Limit ScanRange (nm):	N/A N/A	reak Coverage (%):	N/A



Component Cal Level Table			
Cal Level	Amount		
1	5.000		
2	10.000		
3	25.000		
4	50.000		
5	250.000		
6	500.000		
7	1000.000		
8	1250.000		

Component QC Level Table				
QC Level	Amount			
1	25.000			
2	50.000			
. 3	250.000			
4	500.000			
ICV	600.000			

ICV & CCV Result Table					
Sample ID	Data File Name	Calculated Amount	Area	% Diff	
CALI	A11165002_04	6.43126ug/kg	50181.01	. 28.63	
CAL2	A11165002 05	12.69646ug/kg	112117.31	26.96	
CAL3	A11165002 06	18.30331ug/kg	167545.26	-26.79	
CAL4	A11165002 07	39.95054ug/kg	381544.38	-20.10	
CAL5	A11165002 08	215.03372ug/kg	2112373.87	-13.99	
CAL6	A11165002 09	507.52633ug/kg	5003884.12	1.51	
CAL7	A11165002 10	1028.62353ug/kg	10155323.34	2.86	
CAL8	A11165002 11	1261.43485ug/kg	12456839.09	0.91	
CCV1	A11165002 15	18.54822ug/kg	169966.34	-25.81	
CCV2	A11165002 19	42.04457ug/kg	402245.48	-15.91	
CCV3	A11165002 25	220.25333ug/kg	2163973.65	-11.90	
CCV4	A11165002 26	501.14178ug/kg	4940768.03	0.23	
ICV	A11165002_20	569.27357ug/kg	5614302.18	-5.12	


CANA SEE SEE


Lancaster Laboratories

LCMSMS ANALYSIS REPORT

Component Name:

Hydrazine





Identification Filter:	+ c APCI SRM ms2 209.07@cid20.00 [77.33-77.33, 106.22-106.22]	Component Name: 1st Trace Type:	Hydrazine TIC
2nd Trace Type: Mass Range 2 (m/z): Base Peak(BP):	N/A	Mass Range 1 (m/z): Wavelength Range 2 (nm):	N/A
Retention Time Window (sec): RT Reference: Adjust Using:	30.00000 No N/A	Expected RT (min): View Width (min): Adjust Expected RT:	8.50000 2.00000 No
Detection Options		Peak Detection Algorithm: ICIS Peak Integration	ICIS
ICIS Smoothing Points: Area Noise Factor: ICIS Constrain Peak Width: ICIS Tailing Factor:	5 5 No N/A	Baseline Window: Peak Noise Factor: ICIS Peak Height (%):	100 10 . N/A
ICIS Peak Detection ICIS Minimum Peak Height (S/N): ICIS Window %:	50.0 N/A	ICIS Identify By: ICIS Ion Ratio Confirmation: ICIS Qualifier Ion Coelution (min): ICIS Spectrum Thresholds	Nearest RT N/A N/A
ICIS Forward: ICIS Match:	N/A N/A	ICIS Reverse:	N/A
ICIS Advanced Parameters Minimum Peak Width: Area Tail Extension:	3 5	Noise Method: Multiplet Resolution: Area Scan Window: Calibration	Incos 10 0
Component Type:	Target Compound	%RSD Calculation Method: Internal Standard	Use calculated amounts
ISTD Amount:	N/A	ISTD Units: Target Compounds	N/A
ISTD: Origin: Calibration Curve: Number of Cal, Levels:	IgnoreOrigin Linear 8	Weighting: Response: Target Units: Number of QC Levels: Peak Purity Options	OneOverX Area ug/kg 5
Scan Threshold (mAU): Limit ScanRange (nm):	N/A N/A	Peak Coverage (%):	N/A

Page 5 of 6 Friday, July 01, 2011, 15:15:47



Component Cal Level Table				
Cal Level	Amount			
1	1.000			
2	2.000			
3	5.000			
4	10.000			
5	50.000			
6	100.000			
7	200.000			
8	250.000			

Component QC Level Table				
QC Level	Amount			
1	5.000			
2	10.000			
3	50.000			
4	100.000			
ICV	120.000			

<u>ICV</u>	& C	<u>CCV</u>	Kesul	<u>t lable</u>
Data	File	Nam	e Ca	dculated

Data File Name	Calculated Amount	Агеа	% Diff		
A11165002 04	1.16975ug/kg	16210.29	16.98		
A11165002_05	1.93289ug/kg	30287.34	-3.36		
A11165002 06	4.55090ug/kg	78579.41	8.98		
A11165002 07	9.81926ug/kg	175760.58	-1.81		
A11165002 08	48.77498ug/kg	894343.94	-2.45		
A11165002_09	98.35341ug/kg	1808875.65	-1.65		
A11165002_10	199.06698ug/kg	3666654.19	-0.47		
A11165002_11	254.33183ug/kg	4686078.54	1.73		
A11165002_15	5.29948ug/kg	92387.90	5.99		
A11165002_19	9.26236ug/kg	165487.88	-7.38		
A11165002 25	50.42795ug/kg	924834.87	0.86		
A11165002 26	101.95675ug/kg	1875343.42	1.96		
A11165002_20	121.91427ug/kg	2243482.95	1.60		
	Data File Name A11165002_04 A11165002_05 A11165002_06 A11165002_07 A11165002_09 A11165002_10 A11165002_11 A11165002_15 A11165002_19 A11165002_25 A11165002_26	A11165002_04	Data File Name Calculated Amount Area A11165002_04 1.16975ug/kg 16210.29 A11165002_05 1.93289ug/kg 30287.34 A11165002_06 4.55090ug/kg 78579.41 A11165002_07 9.81926ug/kg 175760.58 A11165002_08 48.77498ug/kg 894343.94 A11165002_09 98.35341ug/kg 1808875.65 A11165002_10 199.06698ug/kg 3666654.19 A11165002_11 254.33183ug/kg 4686078.54 A11165002_15 5.29948ug/kg 92387.90 A11165002_19 9.26236ug/kg 165487.88 A11165002_25 50.42795ug/kg 924834.87 A11165002_26 101.95675ug/kg 1875343.42		

Page 6 of 6 Friday, July 01, 2011, 15:15:47



Sample Name:

CAL1

Data File:

A11165002_04

Std Bracket

Acquisition Date: Sample ID:

06/15/11 06:59:56 PM

Sample Type:

9.98

Vial:

CAL₁ A:3

Run Time(min): Injection Volume(µl): Dilution Factor:

5.00

Instrument Software Version:

Original Data Path:

1.4.1

Instrument Model:

1.00

Instrument Name: Ouantum Instrument Serial Number:

Instrument Method:

TSO Quantum Access C:\XCalibur\Hydrazine TQU01408

C:\XCalibur\Hydrazine

Analysis\Hydraz_soil


Analysis\2011June


Operator:


Quantum


Quan Peak Table


	V	·		, , , , , , , , , , , , , , , , , , ,
Component Name	Calculated Amount	Units	Response Ratio	RT
Monomethylhydrazine	6,409	ug/kg	58896.458	3.71
1,1-Dimethylhydrazine	6.441	ug/kg	50181.009	5.87
Hydrazine	1.170	ug/kg	16210.294	8.43














Page 1 of 1 Thursday, June 16, 2011, 11:37:08



Sample Name:

CAL2

Data File:

A11165002_05

Sample Type:

Std Bracket

Acquisition Date:

06/15/11 07:20:12 PM

Sample ID: Vial:

CAL₂

Run Time(min): Injection Volume(µl):

9.99 5.00

Instrument Software Version:

Original Data Path:

A:4

Dilution Factor:

1.00

Instrument Name:

Instrument Serial Number:

1.4.1

Instrument Model:

TSQ Quantum Access

Ouantum TOU01408

Instrument Method:

C:\XCalibur\Hydrazine

C:\XCalibur\Hydrazine

Analysis\Hydraz_soil


Analysis\2011June


Operator:


Quantum


Ouan Peak Table


		Quant car	· I HAVE		
· —	- Component Name	Calculated Amount	Units	Response Ratio	RT
	Monomethylhydrazine	11.351	ug/kg	105001.171	3.72
	1.1-Dimethylhydrazine	12.699	ug/kg	112117.312	5.89
	Hydrazine	1.933	ug/kg	30287.338	8.47
















Sample Name:

CAL3

Data File:

A11165002_06

Sample Type: Run Time(min):

Std Bracket

5.00

Injection Volume(µl): Dilution Factor:

1.00

Instrument Model: Instrument Method: TSO Quantum Access C:\XCalibur\Hydrazine

Analysis\Hydraz soil

Quantum

Operator:

9.99

Acquisition Date: Sample ID:

06/15/11 07:40:25 PM

CAL₃

A:5

Instrument Software Version: 1.4.1

Instrument Name:

Ouantum

Instrument Serial Number:

TOU01498

Original Data Path:


C:\XCalibur\Hydrazine

Analysis\2011June


Quan Peak Table


Vial:


	Quan 1	CAR LADIC		
- Component Name	Calculated Amount	Units	Response Ratio	RT
Monomethylhydrazine	19.772	ug/kg	183559.444	3.77
1,1-Dimethylhydrazine		ug/kg	167545.255	5.94
Hydrazine	4.551	ug/kg	78579.410	8.52



















Sample Name:

CAL4

Data File:

Sample Type:

A11165002 07

Std Bracket 9.99

Run Time(min): Injection Volume(µl):

5.00

Dilution Factor: Instrument Model:

Instrument Method:

1.00 TSQ Quantum Access C:\XCalibur\Hydrazine

Analysis\Hydraz soil Quantum

Operator:

Acquisition Date:

06/15/11 08:00:42 PM

Sample ID: CAL4

A:6 Vial:

1.4.1 Instrument Software Version:

Instrument Name:

Quantum


Instrument Serial Number: Original Data Path:


TQU01408 C:\XCalibur\Hydrazine


Analysis\2011June


Ouan Peak Table


	VEMI I VIII	* * * * * * * * *		
. Component Name	Calculated Amount	Units -	Response Ratio	RT
 Monomethylhydrazine	44.628	ug/kg	415437.590	3.77
1,1-Dimethylhydrazine	39.925	ug/kg	381544.385	5.92
Hydrazine	9.819	ug/k g	175760.585	8.51














Page 1 of 1 Thursday, June 16, 2011, 11:37:11



Sample Name:

CAL5

Data File:

A11165002 08

Sample Type: Run Time(min): Std Bracket

9.98

5.00

Quantum

Dilution Factor:

Injection Volume(µl):

Instrument Model: Instrument Method: 1.00

TSQ Quantum Access C:\XCalibur\Hydrazine

Analysis\Hydraz soil

Operator:

Acquisition Date:

06/15/11 08:21:00 PM

Sample ID:

CAL5 A:7 1.4.1

Vial: Instrument Software Version:

Instrument Name:

Instrument Serial Number:

Ouantum TQU01408


Original Data Path:

C:\XCalibur\Hydrazine


Analysis\2011June


Onan Peak Table

Vulli 1 cm 1 1 2010					
Component Name	Calculated Amount	Units	Response Ratio	-	RT
Monomethylhydrazine	222.267	ug/kg	2072605.077		3.74
1,1-Dimethylhydrazine	214.824	ug/kg	2112373.869		5.92 -
Hydrazine	48.775	ug/kg	894343.943		8.52





ME_Z/SES Base Peak m/z= 103.50-104.50 F: + c APCI SRM mrz 135.150@dd20.00 [77.325-77.335, 104.135-104.145] MS ICIS A11185002_08





A1185002_06 9470_RT; 5.82_AV; 1_NL; 2.79E5 F; + c APCI SRM mi2 149.100@cd30.00 [77.325-77.335, 106.215-106.225] 105.220 m/z 77,330





Page 1 of 1 Thursday, June 16, 2011, 11:37:12



Sample Name:

CAL₆

Data File:

A11165002 09

Sample Type: Run Time(min): Std Bracket

9.99

Injection Volume(µl): Dilution Factor:

5.00

Instrument Model:

1.00

Instrument Method:

TSQ Quantum Access C:\XCalibur\Hydrazine

Analysis\Hydraz_soil

Quantum

Acquisition Date:

Sample ID:

Vial:

06/15/11 08:41:13 PM

CAL₆

A:8

Instrument Software Version:

1.4.1

Instrument Name:

Quantum

Instrument Serial Number: Original Data Path:

TQU01408


C:\XCalibur\Hydrazine


Analysis\2011June


Operator:


Ouan Peak Table


 Component Name	Calculated Amount	Units	Response Ratio	RT
 Monomethylhydrazine	482.948	ug/kg	4504468.743	3.69
1,1-Dimethylhydrazine	507.010	ug/kg	5003884.123	5.86
Hydrazine	98.353	ug/kg	1808875.647	8.44














Page 1 of I Thursday, June 16, 2011, 11:37:13



Sample Name:

CAL7

Data File:

A11165002_10

Sample Type: Run Time(min):

Std Bracket

9.99

Injection Volume(µI): Dilution Factor:

5.00

Instrument Model: Instrument Method:

1.00

TSQ Quantum Access C:\XCalibur\Hydrazine

Analysis\Hydraz soil

Acquisition Date:

Sample ID:

06/15/11 09:01:30 PM

CAL7

A:9 1.4.1

Instrument Software Version:

Instrument Name:

Instrument Serial Number:

Original Data Path:

Ouantum TOU01408

C:\XCalibur\Hydrazine

Analysis\2011June


Operator:


Quantum


Onan Peak Table


Vial:


V VIII I VIII VIII I VIII I VIII VII					
Component Name	Calculated Amount	Units	Response Ratio	_ RT	
Monomethylhydrazine	1012.909	ug/kg	9448410.120	3.72	
1,1-Dimethylhydrazine	1027.561	ug/kg	10155323.337	5.89	
Hydrazine	199.067	ug/kg -	3666654.195	- 8.4 6	














Page 1 of 1 Thursday, June 16, 2011, 11:37:14



Sample Name:

CAL8

Data File:

A11165002 11

Sample Type: Run Time(min):

Std Bracket

9.99

Injection Volume(µl): Dilution Factor:

5.00 1.00

Instrument Model: Instrument Method: TSQ Quantum Access C:\XCalibur\Hydrazine

Analysis\Hydraz soil

Operator:

Quantum

Acquisition Date:

06/15/11 09:21:47 PM

CAL8

Sample ID:

A:10 1.4.1

Instrument Software Version:

Original Data Path:

Instrument Name:

Quantum

Instrument Serial Number:

TOU01408


C:\XCalibur\Hydrazine


Analysis\2011June


Ouan Peak Table


Vial:


Quali i cak labic:					
Component Name	Calculated Amount	Units	Response Ratio	RT	
Monomethylhydrazine	1289.716	ug/kg	12030709.138	3.72	
1,1-Dimethylhydrazine	1263.240	ug/kg	12487626.542	5.91	
Hydrazine	254.332	ug/kg	4686078.544	8.47	











Page 1 of 1 Thursday, June 16, 2011, 11:37:14



Sample Name:

SYS(MDL)

Data File:

A11165002 03

Sample Type:

Unknown

Run Time(min):

9.99

Injection Volume(µl): Dilution Factor:

5.00

Instrument Model:

1.00

Instrument Method:

TSO Quantum Access C:\XCalibur\Hydrazine

Analysis\Hydraz soil

Operator:

Quantum

Acquisition Date:

06/15/11 06:39:39 PM

SYS(MDL)

Sample ID: Vial:

Instrument Software Version:

A:2 1.4.1

Instrument Name:

Original Data Path:

Quantum

Instrument Serial Number:


TOU01408


C:\XCalibur\Hydrazine


Analysis\2011June


Onan Paak Table


Quali I tak Tabit					
Component Name	Calculated Amount	Units	Response Ratio	RT	
- Monomethylhydrazine	2.772	ug/kg	24970.595	3.74	
1.1-Dimethylhydrazine	3.838	ug/kg	24427.879	5.91	
Hydrazine	0.744	ug/kg	8361.611	8.47	















Page 1 of 1 Thursday, June 16, 2011, 11:37:20



Sample Name:

ICV

Data File:

Operator:

A11165002_20

Sample Type:

Unknown

1.00

Run Time(min):

Dilution Factor:

9.99 Injection Volume(µI): 5.00

Instrument Model: Instrument Method: TSQ Quantum Access

C:\XCalibur\Hydrazine

Analysis\Hydraz soil

Quantum

Acquisition Date:

Sample ID:

06/16/11 12:44:32 AM

ICV

1.4.1

Vial: a:13

Instrument Software Version:

Instrument Name: Instrument Serial Number: Quantum

Original Data Path:


TOU01408


C:\XCalibur\Hydrazine


Analysis\2011June


Quan Peak Table


Component Name	Calculated Amount	Units	Response Ratio	RT
Monomethylhydrazine	522.695	ug/kg	4875265.142	3.67
I,1-Dimethylhydrazine	568.693	ug/kg	5614302.178	5.85
Hydrazine	121.914	ug/kg	2243482.946	8.41














Page 1 of 1 Thursday, June 16, 2011, 11:37:23



Sample Name:

Sample Type:

Run Time(min):

Dilution Factor:

CCV1

Data File:

A11165002 15

OC 9.99

5.00

1.00

Instrument Model: Instrument Method:

Injection Volume(µl):

TSQ Quantum Access

C:\XCalibur\Hydrazine Analysis\Hydraz soil

Acquisition Date: Sample ID:

06/15/11 11:03:08 PM

CCV1

a:5

1.4.1 Instrument Software Version:

Instrument Name:

Instrument Serial Number: Original Data Path:

Quantum TOU01408

C:\XCalibur\Hydrazine

Analysis\2011June


Operator:

Quantum


Quan Peak Table


Vial:


VUMILITAR I ADIO						
Component Name	Calculated Amount	Units	Response Ratio	RT		
Monomethylhydrazine	21.965	ug/kg	204017.632	- 3.71		
1.1-Dimethylhydrazine	18.545	ug/kg	169966.337	5.86		
Hydrazine	5.299	ug/kg	92387.898	8.44		





104 135-104 1451 MS IC25












Page 1 of 1 Thursday, June 16, 2011, 11:37:15



Sample Name:

CCV2

Data File:

A11165002_19.

Acquisition Date:

06/16/11 12:24:15 AM

Sample Type:

QC

Sample ID:

CCV2

Run Time(min):

9.98

Vial: Instrument Software Version: a:6

Injection Volume(µl): Dilution Factor:

5.00 1.00

Instrument Name:

Original Data Path:

1.4.1

Instrument Model:

TSO Quantum Access

Ouantum Instrument Serial Number:

Instrument Method:

C:\XCalibur\Hydrazine

TQU01408

Analysis\Hydraz_soil

C:\XCalibur\Hydrazine


Analysis\2011June


Operator:


Quantum


Ouan Peak Table


·	Component Name	Calculated Amount	Units	Response Ratio		RT
	Monomethylhydrazine	42.810	ug/kg	398473.082	,	3.68
	1,1-Dimethylhydrazine	42.017	ug/kg	402245.481		5.85
	Hydrazine	9.262	ug/kg	165487.879		8.41














Page 1 of 1 Thursday, June 16, 2011, 11:37:16



Sample Name:

CCV3

Data File:

A11165002 25

Acquisition Date:

06/16/11 02:25:56 AM

Sample Type:

QC

Sample ID:

CCV3

Run Time(min):

9.99

Vial:

Original Data Path:

a:7

Injection Volume(µl):

5.00

Instrument Software Version: Instrument Name:

1.4.1

Dilution Factor: Instrument Model: 1.00 TSO Quantum Access

Quantum Instrument Serial Number:

Instrument Method:

C:\XCalibur\Hydrazine

TOU01408

C:\XCalibur\Hydrazine

Analysis\Hydraz soil


Analysis\2011June


Operator:


Quantum


Onan Peak Table


		4 ~	 			
·	Component Name	Calculated Amount	Units	Response Ratio	RT	
	Monomethylhydrazine	222.707	ug/kg	2076710.458	3.74	
	1,1-Dimethylhydrazine	220.038	ug/kg	2163973.647	5.91	
	Hydrazine	50.428	ug/kg	924834.870	8.48	














Page 1 of 1 Thursday, June 16, 2011, 11:37:17

Raw QC Data



Sample Name:

BLK Sand

Data File:

A11165002_14

Sample Type:

Blank

Run Time(min): Injection Volume(µl): 9.99 5.00

Dilution Factor: Instrument Model:

1.00

Instrument Method:

TSQ Quantum Access C:\XCalibur\Hydrazine

Analysis\Hydraz_soil

Operator:

Acquisition Date: Sample ID:

06/15/11 10:22:36 PM

BLK. Sand

Vial:

Instrument Software Version:

1.4.1

a:11

Instrument Name:

Instrument Serial Number:

Ouantum TQU01408

Original Data Path:


C:\XCalibur\Hydrazine


Analysis\2011June


Quantum


Quan Peak Table


V 4411 1 4011 1 4011 1 4011 1 4011 1 4011 1 4011 1 4011 1 4011 1 4011 1 4011 1 4011 1 4011 1 4011 1 4011 1 4011					
	Component Name	Calculated Amount	Units	Response Ratio	RT
	Hydrazine	N/A	ug/kg	N/A	- N/A
	1,1-Dimethylhydrazine	N/A	ug/kg	N/A	,N/A
	Monomethylhydrazine	. N/A	ug/kg	N/A	N/A














Page 1 of 1 Thursday, June 16, 2011, 11:37:19



Sample Name:

6310730 (MS)

Data File:

Operator:

A11165002 23

Sample Type:

Unknown

Run Time(min):

9.99 5.00

Injection Volume(µI):

Dilution Factor:

Instrument Model:

Instrument Method:

1.00

TSQ Quantum Access C:\XCalibur\Hydrazine

Analysis\Hydraz_soil

Quantum

Acquisition Date:

Sample ID:

Vial: Instrument Software Version:

Instrument Name:

Instrument Serial Number:

Original Data Path:

06/16/11 01:45:22 AM

6310730 (MS)

a:16

1.4.1

Quantum


TOU01408


C:\XCalibur\Hydrazine


Analysis\2011June


Ouan Peak Table


V UNIT VALLE I AND					
Component Name	Calculated Amount	Units	Response Ratio	.	RT
Monomethylhydrazine	50.613	ug/kg	471272.236		3.69
1,1-Dimethylhydrazine	55.446	ug/kg	535146.490	:	5.86
Hydrazine	13.766	ug/kg	248567.911	-	8.43














Page 1 of 1 Thursday, June 16, 2011, 11:37:25



Sample Name:

6310731(MSD)

Data File:

A11165002 24

Sample Type:

Unknown

Run Time(min):

9.99

Injection Volume(µl): Dilution Factor:

5.00

Instrument Model: Instrument Method: 1.00

TSQ Quantum Access C:\XCalibur\Hydrazine

Analysis\Hydraz soil

Acquisition Date:

Sample ID:

06/16/11 02:05:39 AM

6310731(MSD)

Instrument Software Version:

Instrument Name:

Quantum

a:17

1.4.1

Instrument Serial Number: Original Data Path:

TOU01408

C:\XCalibur\Hydrazine


Analysis\2011June


Quantum Operator:


Ouan Peak Table


Vial:


 Component Name	Calculated Amount	Units	Response Ratio	RT
 Monomethylhydrazine	36.384	ug/kg	338527.471	3.71
1,1-Dimethylhydrazine	55.836	ug/kg	539001.853	5.90
Hydrazine	12.563	ug/kg	226365.897	8.46














Page 1 of 1 Thursday, June 16, 2011, 11:37:26



Sample Name:

LCS

Data File:

Sample Type:

Run Time(min): Injection Volume(µI):

Dilution Factor:

Instrument Model:

Instrument Method:

A11165002 21

Unknown 9.98

5.00

1.00 TSQ Quantum Access

C:\XCalibur\Hydrazine Analysis\Hydraz_soil

Quantum

Acquisition Date:

Sample ID:

06/16/11 01:04:49 AM

LCS a:14

Instrument Software Version: 1.4.1

Instrument Name:

Original Data Path:

Instrument Serial Number:

Ouantum

TQU01408

C:\XCalibur\Hydrazine


Analysis\2011June


Operator:


Ouan Peak Table


Vial:


Component Name	Calculated Amount	Units	Response Ratio	RT
Monomethylhydrazine	424.297	ug/kg	3957319.268	3.69
1,1-Dimethylhydrazine	516.922	ug/kg	5101969.819	5.88
Hydrazine	92.056	ug/kg	1692721.302	8.45














Page I of 1 Thursday, June 16, 2011, 12:53:34



Sample Name:

LCSD

Data File:

A11165002 22

Sample Type:

Unknown

Run Time(min): Injection Volume(µl): 9.98 5.00

Dilution Factor:

1.00

Instrument Model: Instrument Method:

TSQ Quantum Access C:\XCalibur\Hydrazine

Analysis\Hydraz_soil

Quantum Operator:

Acquisition Date:

06/16/11 01:25:06 AM

LCSD

Sample ID:

a:15 1.4.1

Instrument Software Version:

Instrument Name:

Original Data Path:

Instrument Serial Number:

Quantum TQU01408


C:\XCalibur\Hydrazine


Analysis\2011June


Ouan Peak Table


Vial:


Component Name	Calculated Amount	Units	Response Ratio	RT
Monomethylhydrazine	481.020	ug/kg	4486479.012	3.71
1,1-Dimethylhydrazine	548.695	ug/kg	5416398.549	-i- 😅 5.88
Hydrazine	113.592	ug/kg	2089963.356	- 8.45














Page 1 of 1 Thursday, June 16, 2011, 11:37:25

Preparation Logs

Organic Extraction Batchlog Assigned to: 2628 Meng Yu

Reviewed by: Of mo-

Start Date: 6 15111 Tech 1: Any 2.6 2.8

Tech 2:

Start time: _

ر اران اران

11165002

Comments 34°C ł hA-7 <u>జ</u> 핊 Hydrazines in Soil 표 Amt FV (mL.) (mL.) 100 <u>ئ</u> -ار د <u>ط</u> = ت کا ا ان o, O **↓**2 12882-4-1 0.10 0.10 125952-19-1 MS Sol. Ϋ́Х Amt (mL) **3** با مح 72 Ł 12 SS/IS Sol. **∤** ب<u>ا</u> 2 ₹ **∤** 2 بد 2 Sample Amt Code (q) Prep Analysis: 00000 OPR165002 BLK165002 448-1 448-1 6310731MSD 6310730MS Dept: 37 BLANKA LCSDA မွ LCSA

Prio	٥	۲.	Ъ		
Due Date	1 100 00 00	100/21/201	06/21/2011		
Analyses		10346	10346		
Comments					
BC	1	34°a	र रु		
듄					
Ŧ		'n	Ś		
₹ (E	1	12 12	15.1 XV		
Amt (m)	Ì	1 2	苕		
SS/IS Sal.		ŁŻ	松		
Amt	2		_		
Sample	COOE	448-D	448-1		
Sample #		1 6310728	2 6310729BKG		

0L#72

Work Station	S-bath ID		C S-bath ID
Balance #	Documented	temps are	Documented temps are NIST corrected.

11165002

ပ

C M-vap

C N-Evap

Page 1 of 1

DF = Dilution Factor FV = Final Volume

Internal Standard

Rack ID:

Moisture Data



CLIENT: Olin Corporation

SDG: OLN72

SAMPLE NUMBERS:

Sample #	Sample Code
6310728	448-DFD
6310729	448-1BKG
6310730	448-1MS
6310731	448-1MSD

Laboratory Compliance Quality Control

Analysis Name	LCS <u>%REC</u>	LCSD <u>%REC</u>	LCS/LCSD <u>Limits</u>	<u>RPD</u>	RPD Max
Batch number: 11165820002B Moisture Moisture Moisture Duplicate	Sample 100 100 100	number(s): 6310728 99-101 99-101 99-101	- 6 31073	1

Sample Matrix Quality Control

Analysis Name	BKG Conc	DUP <u>Conc</u>	RPD	RPD Max
Batch number: 11165820002B	Sample	e number(s): 6310	728-6310731
Moisture	6.5	6.1	7	15
Moisture Duplicate	6.5	6.1	7	15
Moisture	6.5	6.1	7	15

^{* -} Outside of specification

Moisture Data Report

Batch #: 11165820002

				Sample			Analysis	verified
Sample ID	Batch ID	Analysis#	Tare Wt	<u>Wt</u>	Dry Wt	<u>%Moisture</u>	Date (Emp#)	Date (Emp#)
LCS 89.5% Std.			1.0844	5.0155	1.6177	89.37) 6/15/11 (1382/SAS)
6310728FD	В	00111	1.0543	9.3679	9.8431	6.18) 6/15/11 (1382/SAS)
6310729BKG	В	00111	1.0750	9.2926	9.7591	6.55	•) 6/15/11 (1382/SAS)
6310730MS	В	00118				6.55	•) 6/15/11 (1382/SAS)
6310731MSD	В	00118				6.55	•) 6/15/11 (1382/SAS)
6310731MSD	В	00121	1.0689	9.2990	9.7989	6.12	6/14/11 (1201/SWF) 6/15/11 (1382/SAS)

OLHVZ BRAG

^{(1) -} The result for one or both determinations was less than five times the LOQ.