

PROHEXADIONE-CALCIUM (BX-112)

014083

Ames assay (§84-2)

EPA Reviewer: Irving Mauer, PhD
Registration Action Branch 2/HED (7509C)
EPA Work Assignment Manager: Marion Copley, DVM
Toxicology Branch 1/HED (7509C)

Irving Mauer 8/27/99
Marion Copley 8/27/99

DATA EVALUATION RECORD

STUDY TYPE: Bacterial reverse mutation test

OPPTS Number: 870.5100

OPP Guideline No.: §84-2

DP BARCODE: D246707

SUBMISSION CODE: S543930

P.C. CODE: 112600

TOX. CHEM. NO.: None

TEST MATERIAL (PURITY): BX-112 Technical (89.8% a.i.)

SYNONYMS: Calcium salt of 3,5-dioxo-4-propionyl-cyclohexane-1-carboxylic acid

CITATION: Jones, E. and L.A. Wilson. (1997) Bacterial Reverse Gene Mutation Assay of BX-112 Technical. Huntingdon Research Centre Ltd., Huntingdon, Cambridgeshire, PE 18 6ES, England. HRC Project No. KCI 25/881559, BASF Registration Document No. 91/10807. MRID 44457765. Unpublished.

SPONSOR: BASF Corporation, Agricultural Products, P.O. Box 13528, Research Triangle Park, NC

EXECUTIVE SUMMARY:

In two trials of a reverse gene mutation assay in bacteria (MRID 44457765), *S. typhimurium* strains TA 98, TA 100, TA 1535, TA 1537, TA 1538, and *E. coli* strain WP2 uvrA were exposed to prohexadione-calcium (BX-112 Technical, 89.8% a.i.) in water, at concentrations of 312.5, 625, 1,250, 2,500, and 5,000 µg/plate in the presence and absence of mammalian metabolic activation. Appropriate solvent and positive controls were included.

Prohexadione-calcium was tested up to the limit concentration, 5,000 µg/plate, and was not toxic to strains TA 98, TA 100, TA 1535, TA 1537, TA 1538, or WP2 uvrA in the absence or presence of metabolic activation. Prohexadione-calcium was not mutagenic to *S. typhimurium* strains TA 98, TA 100, TA 1535, TA 1537, TA 1538, or WP2 uvrA with or without metabolic activation at up to 5,000 µg/plate. The positive controls did induce the appropriate responses in the corresponding strains.

This study is classified as **acceptable (§84-2)**. It does satisfy the requirement for FIFRA Test Guideline for in vitro mutagenicity (bacterial reverse gene mutation) data.

PROHEXADIONE-CALCIUM (BX-112)

Ames assay (§84-2)

COMPLIANCE: Signed and dated GLP, Quality Assurance, and Data Confidentiality statements were provided.

I. MATERIALS AND METHODS

A. MATERIALS

1. Test Material: BX-112 Technical

Description: Pale yellow powder


Lot/Batch #: G14-04

Purity: 89.8% a.i.

Stability of compound: Expiration date reported to be September 1989; stored at 4°C in the dark. The study was conducted between September 20, 1988 and October 10, 1988.

CAS #: 127277-53-6

Structure:



Solvent used: Water

Other comments: The test substance did not dissolve in any of the solvents tested but formed a stable suspension in water.

2. Control Materials:

Negative: Solvent served as negative control

Solvent/final concentration: Water/0.1 mL/plate

Positive: Nonactivation:

Sodium azide NA^a $\mu\text{g/plate}$ TA100, TA15352-Nitrofluorene 1 $\mu\text{g/plate}$ TA982-Nitrofluorene 2 $\mu\text{g/plate}$ TA15389-Aminoacridine 80 $\mu\text{g/plate}$ TA1537

Other: NA

N-ethyl-N'-nitro-N-nitrosoguanidine 3 $\mu\text{g/plate}$ TA100N-ethyl-N'-nitro-N-nitrosoguanidine 2 $\mu\text{g/plate}$ WP2 uvrAN-ethyl-N'-nitro-N-nitrosoguanidine 5 $\mu\text{g/plate}$ TA 1535^a NA = Not Applicable

Activation:

2-Aminoanthracene (2-anthramine) 0.5 $\mu\text{g/plate}$ *S. typhimurium* TA98 and TA15382-Aminoanthracene 1 $\mu\text{g/plate}$ *S. typhimurium* TA1002-Aminoanthracene 2 $\mu\text{g/plate}$ *S. typhimurium* TA1535 and TA1537

2-Aminoanthracene 20 $\mu\text{g}/\text{plate}$ *E. coli* WP2 uvrA

Other (list): NA

3. Activation: S9 derived from

<input checked="" type="checkbox"/> Aroclor 1254	<input checked="" type="checkbox"/> induced	<input checked="" type="checkbox"/> rat	<input checked="" type="checkbox"/> liver
<input type="checkbox"/> phenobarbital	<input type="checkbox"/> non-induced	<input type="checkbox"/> mouse	<input type="checkbox"/> lung
<input type="checkbox"/> none	<input type="checkbox"/> hamster	<input type="checkbox"/> other	
<input type="checkbox"/> other	<input type="checkbox"/> other		

The S9 mix was derived from rats pretreated with Aroclor 1254 and 1 mL of S9 mix contained the following: 0.1 mL S9 fraction (10% v:v), 8 μmol MgCl_2 , 33 μmol KCl, 4 μmol NADH, 4 μmol NADPH, 5 μmol glucose-6-phosphate, and 100 μmol phosphate buffer (pH 7.4). The final concentration of S9 in culture was approximately 2%.

4. Test organisms: *S. typhimurium* strains☐ TA97 ☒ TA98 ☒ TA100 ☐ TA102 ☐ TA104☒ TA1535 ☒ TA1537 ☒ TA1538 ; list any others: *E. coli* WP2 uvrA

Properly maintained? Yes

Checked for appropriate genetic markers (rfa mutation, R factor)? Yes

5. Test compound concentrations used (preliminary cytotoxicity test, if performed and main assay):a. Preliminary cytotoxicity test

Nonactivated conditions: 5, 50, 500, and 5,000 $\mu\text{g}/\text{plate}$ TA 98, TA 100, TA 1535, TA 1537, TA 1538, WP2 uvrA, one plate/dose

Activated conditions: 5, 50, 500, and 5,000 $\mu\text{g}/\text{plate}$ TA 98, TA 100, TA 1535, TA 1537, TA 1538, WP2 uvrA, one plate/dose

b. Gene mutation assay

Nonactivated conditions: 312.5, 625, 1,250, 2,500, and 5,000 $\mu\text{g}/\text{plate}$ TA 98, TA 100, TA 1535, TA 1537, TA 1538, WP2 uvrA, three plates/dose

Activated conditions: 312.5, 625, 1,250, 2,500, and 5,000 $\mu\text{g}/\text{plate}$ TA 98, TA 100, TA 1535, TA 1537, TA 1538, WP2 uvrA, three plates/dose

B. TEST PERFORMANCE1. Type of Salmonella assay:☒ standard plate test

237

- pre-incubation (— minutes)
- "Prival" modification (i.e. azo-reduction method)
- spot test
- other (describe)

2. Protocol

a. Preliminary cytotoxicity test

A 0.1 mL aliquot of a bacterial culture containing approximately 2×10^9 cells/mL and 0.5 mL S-9 mix or 0.5 mL 0.1 M sodium orthophosphate buffer were placed in a glass bijou bottle. A 0.1 mL aliquot of the test solution was added to the bottle, followed by 2 mL of histidine or tryptophan deficient agar. The mixture was thoroughly shaken and overlaid onto previously prepared plates containing 25 mL of minimal agar. Single plates were used for each dose level and solvent controls were included. The plates were incubated at 37°C for 72 hours, then examined for the appearance of a complete bacterial lawn and scored for revertant colonies. Toxicity was assessed as a substantial reduction in revertant colony counts or by the absence of a complete background bacterial lawn.

b. Gene mutation assay

A 0.1 mL aliquot of bacterial suspension and 0.5 mL of sterile 0.1 M sodium orthophosphate buffer (pH 7.4) were added to each of one set of sterile bijou bottles (3 bottles/dose). A 0.1 mL aliquot of the test substance was added to each bottle at one of the five test concentrations separated by a two-fold interval. A 2 mL aliquot of histidine or tryptophan deficient agar was added to each bottle, thoroughly mixed, and overlaid onto previously prepared plates containing 25 mL of minimal agar. Three plates were used for each test concentration. Solvent and positive controls were included. The plates were incubated at 37°C for 72 hours. For metabolic activation cultures, 0.5 mL of S9 mix was substituted for the phosphate buffer. Solvent and positive controls were included. The plates were scored for revertant colonies. A positive test resulted if there was a reproducible, statistically significant (not defined), dose-related increase in the number of revertant colonies.

II. REPORTED RESULTS

A. Analytical Determinations

Analytical determinations of the dose formulations were not performed.

B. Preliminary cytotoxicity assay

One trial was conducted with prohexadione-calcium at four concentrations ranging from 5 to 5,000 $\mu\text{g}/\text{plate}$ with and without metabolic activation. The tester strains used were TA 98, TA 100, TA 1535, TA 1537, TA 1538, and WP2 uvrA. Single plates were used per dose. The test compound was not toxic to any of the tester strains. The results of the preliminary cytotoxicity assay are presented in study report Table 1, page 9.

208

C. Mutagenicity assay

Two trials were conducted with prohexadione-calcium at five concentrations ranging from 312.5 to 5,000 $\mu\text{g}/\text{plate}$ with or without metabolic activation (Aroclor1254-induced rat livers). The tester strains used were TA 98, TA 100, TA 1535, TA 1537, TA 1538, and WP2 uvrA. Three plates were used per dose. No mutagenic activity was observed in the solvent control. No substantial increase in mutagenic activity was observed in revertant colonies of any of the six tester strains treated with the test compound at any dose level either with or without metabolic activation. The positive controls gave the appropriate responses. The individual plate results of the gene mutation assay are presented in study report Tables 3, 4, 6, and 7, pages 11, 12, 13, 15, 16, and 17). The results are summarized in study report Tables 2 and 5, pages 10 and 14, and are presented as an Appendix to this DER.

III. REVIEWER'S DISCUSSION/CONCLUSIONS:

- A. Prohexadione-calcium was assayed in six *S. typhimurium* strains with and without S9 at concentrations up to the limit dose. The test compound was not cytotoxic, and failed to induce a genotoxic response. The sensitivity of the test system to detect mutagenesis was adequately demonstrated by the response obtained with the nonactivated and S9-activated positive controls.
- B. Study deficiencies - Data on the analysis of dose formulations for actual concentration were not submitted. However, as the test compound was tested up to the limit dose, this deficiency is not expected to affect the interpretation of this study. Several minor deficiencies noted in the study that are not considered to affect the validity of the study results are:
- Criteria for statistical significance were not provided.
 - Historical control ranges (solvent and positive) were not provided.

DER Bacterial Reverse Mutation

MRIID 44457765

Page is not included in this copy.

Pages 7 through 9 are not included.

The material not included contains the following type of information:

- ☐ Identity of product inert ingredients.
- ☐ Identity of product impurities.
- ☐ Description of the product manufacturing process.
- ☐ Description of quality control procedures.
- ☐ Identity of the source of product ingredients.
- ☐ Sales or other commercial/financial information.
- ☐ A draft product label.
- ☐ The product confidential statement of formula.
- ☐ Information about a pending registration action.
- ☒ FIFRA registration data.
- ☐ The document is a duplicate of page(s)
- ☐ The document is not responsive to the request.

The information not included is generally considered confidential by product registrants. If you have any questions, please contact the individual who prepared the response to your request.