




Agrochemical Analysis by Photoionization MS

**Matthew D. Evans, Karl A. Hanold, and
Jack A. Syage**

FPRW
July 28-31, 2002



Syagen Technology, Inc.
1411 Warner Avenue
Tustin, CA 92780
www.syagen.com



Outline

- Problem
 - Agrochemicals including pesticides are not routinely detectable by ESI and APCI*
 - Chemical weapons have similar detection problems*
 - This work evaluates the applicability of photoionization (PI) MS*
- Experimental
 - Low pressure photoionization (LPPI) with quadrupole ion trap TOF (QitTof™) MS — (Syagen Radiance Pro)*
 - APPI and APCI with flow injection MS — (Agilent 1100 Series LC/MSD)*
- Results
 - Comparison of LPPI to atmospheric pressure sources APPI and APCI*
 - Detectability of CW agents and byproducts in soil extracts and contaminated water*
 - Lessons learned*
- Conclusions



Chlorinated Pesticides: Non-aromatic vs. Aromatic

IP's of non-aromatic chlorinated cyclic hydrocarbons

Cyclohexane: IP = 9.88 eV

Chlorocyclohexane: IP = 10.10 eV




Hexachlorocyclohexane (α -BHC): not observed

Dibutyl Chlorendate: $10 < IP < 10.6$???


IP's of highly chlorinated aromatics

Benzene: IP = 9.24 eV

Hexachlorobenzene: 9.0 eV



LPPI vs. APPI vs. APCI



Characteristics of compounds that photoionize

Aromatics work well

LPPI has higher detection rate than APPI, which in turn is moderately better than APCI

Other observations

Ion dissociation efficient for compounds containing $-CCl_3$ (e.g. p,p'-DDT and methoxychlor)

MH⁺ detected for compounds with basic functional groups

