

METHOD CONTROL PARAMETERS

1

Method Information For: C:\MSDCHEM\1\METHODS\ (b) (7) RP1581_051914.M
Method Sections To Run:

() Save Copy of Method With Data
() Instrument Control Pre-Run Cmd/Macro =
() Data Analysis Pre-Run Cmd/Macro =
(X) Data Acquisition
(X) Data Analysis
() Instrument Control Post-Run Cmd/Macro =
() Data Analysis Post-Run Cmd/Macro =
Method Comments:

END OF METHOD CONTROL PARAMETERS

Quantitation method

(b) (7)(C) 8/11/14

INSTRUMENT CONTROL PARAMETERS: Enhanced

C:\MSDCHEM\1\METHODS\ (b) (7) RP1581_051914.M
Fri Aug 08 13:44:40 2014

Control Information

Sample Inlet : GC
Injection Source : GC ALS
Mass Spectrometer : Enabled

Oven
Equilibration Time 0.5 min
Oven Program On
35 °C for 2 min
then 20 °C/min to 325 °C for 3 min
Run Time 19.5 min

Front Injector
Syringe Size 10 µL
Injection Volume 1 µL
Injection Repetitions 1
Solvent A Washes (PreInj) 2
Solvent A Washes (PostInj) 3
Solvent A Volume 8 µL
Solvent B Washes (PreInj) 2
Solvent B Washes (PostInj) 3
Solvent B Volume 8 µL
Sample Washes 1
Sample Wash Volume 8 µL
Sample Pumps 3
Dwell Time (PreInj) 0 min
Dwell Time (PostInj) 0 min
Solvent Wash Draw Speed 300 µL/min
Solvent Wash Dispense Speed 6000 µL/min
Sample Wash Draw Speed 300 µL/min
Sample Wash Dispense Speed 6000 µL/min
Injection Dispense Speed 6000 µL/min
Viscosity Delay 0 sec
Sample Depth Disabled

Back Injector
Syringe Size 10 µL
Injection Volume 1 µL
Injection Repetitions 1
Solvent A Washes (PreInj) 0
Solvent A Washes (PostInj) 0
Solvent A Volume 8 µL
Solvent B Washes (PreInj) 0
Solvent B Washes (PostInj) 0
Solvent B Volume 8 µL
Sample Washes 0
Sample Wash Volume 8 µL
Sample Pumps 6
Dwell Time (PreInj) 0 min
Dwell Time (PostInj) 0 min
Solvent Wash Draw Speed 300 µL/min
Solvent Wash Dispense Speed 6000 µL/min
Sample Wash Draw Speed 300 µL/min
Sample Wash Dispense Speed 6000 µL/min
Injection Dispense Speed 6000 µL/min
Viscosity Delay 0 sec
Sample Depth Disabled

Front SS Inlet He

Mode	Split	
Heater	On	275 °C
Pressure	On	16.494 psi
Total Flow	On	114.93 mL/min
Septum Purge Flow	On	3 mL/min
Gas Saver	On	30 mL/min After 2 min
Split Ratio	20 :1	
Split Flow	106.6 mL/min	

Back SS Inlet He

Mode	Split	
Heater	Off	
Pressure	Off	
Total Flow	Off	
Septum Purge Flow	Off	
Gas Saver	Off	
Split Ratio	100 :1	
Split Flow	277.4 mL/min	

Column #1

450 °C: 25 m x 320 µm x 0 µm
 In: Front SS Inlet He
 Out: Front Detector TCD

(Initial)	35 °C
Pressure	16.494 psi
Flow	5.3298 mL/min
Average Velocity	70.102 cm/sec
Holdup Time	0.59437 min
Flow Program	Off
5.3298 mL/min for 0 min	
Run Time	19.5 min

Column #2

450 °C: 25 m x 320 µm x 0 µm
 In: Back SS Inlet He
 Out: Front Detector TCD

(Initial)	35 °C
Pressure	10 psi
Flow	2.774 mL/min
Average Velocity	43.399 cm/sec
Holdup Time	0.96009 min
Pressure Program	Off
10 psi for 0 min	
Run Time	19.5 min

Front Detector TCD

Heater	On	300 °C
Reference Flow	Off	
Makeup Flow	On	20 mL/min
Const Col + Makeup	Off	
Negative Polarity	Off	
Filament	Off	

1
Signals
Front Signal Save On
Test Plot Save On
Test Plot Save On
Test Plot Save On

MS ACQUISITION PARAMETERS

General Information

Tune File : dftppnew051614.u
Acquisition Mode : Scan

MS Information

Solvent Delay : 3.50 min

EMV Mode : Relative
Relative Voltage : -71

[Scan Parameters]

Low Mass : 35.0
High Mass : 500.0
Threshold : 50
Sample # : 3 A/D Samples 8
Plot 2 low mass : 50.0
Plot 2 high mass : 550.0

[MSZones]

MS Source : 230 C maximum 250 C
MS Quad : 150 C maximum 200 C

END OF MS ACQUISITION PARAMETERS

TUNE PARAMETERS for SN: <offline>

Trace Ion Detection is OFF.

EMISSION : 34.600
ENERGY : 70.000
REPELLER : 29.960
IONFOCUS : 90.200
ENTRANCE_LE : 0.000
EMVOLTS : 1200.000

Actual EMV : -1
GAIN FACTOR : <Unable to calculate gain facto

r.>

AMUGAIN : 4095.000
AMUOFFSET : 255.000
FILAMENT : 1.000
DCPOLARITY : 0.000
ENTLENSOFFS : 25.100
MASSGAIN : -2048.000
MASSOFFSET : -499.000

END OF TUNE PARAMETERS

END OF INSTRUMENT CONTROL PARAMETERS

DATA ANALYSIS PARAMETERS

Method Name: C:\MSDCHEM\1\METHODS\ (b) (7) (C) RP1581_051914.M

Percent Report Settings

Sort By: Signal

Output Destination

Screen: No

Printer: No

File: No

Integration Events: Meth Default

Generate Report During Run Method: Yes

Signal Correlation Window: 0.020

Qualitative Report Settings

Peak Location of Unknown: Apex

Library to Search Minimum Quality

C:\Database\NIST05a.L 0

Integration Events: Meth Default

Report Type: Summary

Output Destination

Screen: No

Printer: No

File: No

Generate Report During Run Method: Yes

Quantitative Report Settings

Report Type: Detailed

Output Destination

Screen: No

Printer: Yes

File: epatemp.txt

Generate Report During Run Method: No

Freedom Industries

Calibration Last Updated: Wed Aug 06 11:32:39 2014

Reference Window: 2.00 Minutes

Non-Reference Window: 1.00 Minutes

(b) (7) (C) RP1581_051914.M Fri Aug 08 13:44:39 2014

Correlation Window: 0.10 minutes
Default Multiplier: 1.00
Default Sample Concentration: 0.00

Compound Information

1) Cyclohexanol

()

Ret. Time 4.613 min., Extract & Integrate from 4.113 to 5.113 min.

Signal	Rel Resp.	Pct. Unc.(rel)	Integration
Tgt 57.10			*** METH DEFAULT ***
Q1 82.10	57.10	20.0	*** METH DEFAULT ***
Q2 67.10	43.80	20.0	*** METH DEFAULT ***

Lvl ID	Conc (ug/mL)	Response
1	50.000	156539
2	50.000	155813
3	50.000	169255
4	50.000	157595
5	50.000	151352
6	50.000	168333
7	50.000	167964
8	50.000	163415

Qualifier Peak Analysis ON
Curve Fit: Avg. RF

2) MCHM (2 peaks combined)

()

Ret. Time 6.530 min., Extract & Integrate from 6.030 to 7.030 min.

Signal	Rel Resp.	Pct. Unc.(rel)	Integration
Tgt 55.10			*** METH DEFAULT ***
Q1 97.10	39.40	20.0	*** METH DEFAULT ***
Q2 81.10	30.30	20.0	*** METH DEFAULT ***

Lvl ID	Conc (ug/mL)	Response
1	0.470	874
2	0.946	2156
3	2.840	7147
4	4.730	11471
5	9.460	22258
6	47.300	166406
7	94.600	367991
8	189.000	789691

Qualifier Peak Analysis ON
Curve Fit: Linear

3) PPH

()

Ret. Time 7.784 min., Extract & Integrate from 7.284 to 8.284 min.

Signal	Rel Resp.	Pct. Unc.(rel)	Integration
Tgt 94.10			*** METH DEFAULT ***
Q1 77.10	23.70	20.0	*** METH DEFAULT ***
Q2 152.10	18.60	20.0	*** METH DEFAULT ***

Lvl ID	Conc (ug/mL)	Response
1	0.470	604
2	0.946	1506
3	2.840	6257
4	4.730	10474
5	9.460	21573
6	47.300	189499
7	94.600	441478
8	189.000	1023617

1
Qualifier Peak Analysis ON
Curve Fit: Linear

4) 1,4-CHDM

()

Ret. Time 8.592 min., Extract & Integrate from 8.092 to 9.092 min.

Signal	Rel Resp.	Pct. Unc.(rel)	Integration
Tgt 95.10			*** METH DEFAULT ***
Q1 67.10	41.40	20.0	*** METH DEFAULT ***
Q2 41.10	53.20	20.0	*** METH DEFAULT ***

Lvl ID	Conc (ug/mL)	Response
1	0.470	280
2	0.946	545
3	2.840	2435
4	4.730	4360
5	9.460	8585
6	47.300	74923
7	94.600	188094
8	189.000	411739

Qualifier Peak Analysis ON
Curve Fit: Linear

5) DMCH-1,4-DC

()

Ret. Time 9.152 min., Extract & Integrate from 8.652 to 9.652 min.

Signal	Rel Resp.	Pct. Unc.(rel)	Integration
Tgt 81.10			*** METH DEFAULT ***
Q1 108.10	36.20	20.0	*** METH DEFAULT ***
Q2 140.10	63.60	20.0	*** METH DEFAULT ***

Lvl ID	Conc (ug/mL)	Response
1	0.470	646
2	0.946	1398
3	2.840	5114
4	4.730	7748
5	9.460	18343
6	47.300	130770
7	94.600	293609
8	189.000	639126

Qualifier Peak Analysis ON
Curve Fit: Linear

6) di-PPH

()

Ret. Time 9.856 min., Extract & Integrate from 9.356 to 10.356 min.

Signal	Rel Resp.	Pct. Unc.(rel)	Integration
Tgt 59.10			*** METH DEFAULT ***
Q1 94.10	30.00	20.0	*** METH DEFAULT ***
Q2 210.20	9.00	20.0	*** METH DEFAULT ***

Lvl ID	Conc (ug/mL)	Response
1	0.470	401
2	0.946	1017
3	2.840	2992
4	4.730	4799
5	9.460	11129
6	47.300	83676
7	94.600	214860
8	189.000	-1

Qualifier Peak Analysis ON
Curve Fit: Linear

END OF DATA ANALYSIS PARAMETERS

Fri Aug 08 13:44:41 2014

Quantitation Report (QT Reviewed)

Data Path : C:\msdchem\1\DATA\ (b) (7) \RP1581\
Data File : 051914_001.D
Acq On : 19 May 2014 11:42 am
Operator : (b) (7)
Sample : DFTPP
Misc : DFTPP
ALS Vial : 1 Sample Multiplier: 1

Quant Time: Aug 06 11:31:39 2014
Quant Method : C:\msdchem\1\METHODS\ (b) (7) RP1581_051914.M
Quant Title : Freedom Industries
QLast Update : Wed Aug 06 11:28:36 2014
Response via : Initial Calibration

Compound	R.T.	QIon	Response	Conc	Units	Dev (Min)


System Monitoring Compounds						
1) Cyclohexanol	0.000	57	0	0.00	ug/mL	
Target Compounds						
2) MCHM (2 peaks combined)	0.000		0	N.D.		Qvalue
3) PPH	0.000		0	N.D.		
4) 1,4-CHDM	0.000		0	N.D.		
5) DMCH-1,4-DC	0.000		0	N.D.		
6) di-PPH	0.000		0	N.D.		

(#) = qualifier out of range (m) = manual integration (+) = signals summed

Quantitation Report (QT Reviewed)

Data Path : C:\msdchem\1\DATA(b) (7)(C) RP1581\
 Data File : 051914_001.D
 Acq On : 19 May 2014 11:42 am
 Operator : (b) (7)(C)
 Sample : DFTPP
 Misc : DFTPP
 ALS Vial : 1 Sample Multiplier: 1

Quant Time: Aug 06 11:31:39 2014
 Quant Method : C:\msdchem\1\METHODS(b) (7)(C) P1581_051914.M
 Quant Title : Freedom Industries
 QLast Update : Wed Aug 06 11:28:36 2014
 Response via : Initial Calibration



Quantitation Report (QT Reviewed)

Data Path : C:\msdchem\1\DATA (b) (7)(C) RP1581\
Data File : 051914_002.D
Acq On : 19 May 2014 12:17 pm
Operator : (b) (7)(C)
Sample : Blank
Misc : Blank
ALS Vial : 2 Sample Multiplier: 1

Quant Time: Aug 06 11:32:06 2014 (b) (7)
Quant Method : C:\msdchem\1\METHODS (C) RP1581_051914.M
Quant Title : Freedom Industries
QLast Update : Wed Aug 06 11:28:36 2014
Response via : Initial Calibration

Compound	R.T.	QIon	Response	Conc	Units	Dev (Min)


System Monitoring Compounds						
1) Cyclohexanol	0.000	57	0	0.00	ug/mL	
Target Compounds						Qvalue
2) MCHM (2 peaks combined)	0.000		0	N.D.		
3) PPH	0.000		0	N.D.		
4) 1,4-CHDM	0.000		0	N.D.		
5) DMCH-1,4-DC	0.000		0	N.D.		
6) di-PPH	0.000		0	N.D.		

(#) = qualifier out of range (m) = manual integration (+) = signals summed

Quantitation Report (QT Reviewed)

Data Path : C:\msdchem\1\DATA (b) (7)(C) RP1581\
 Data File : 051914_002.D
 Acq On : 19 May 2014 12:17 pm
 Operator : (b) (7)
 Sample : (C) k
 Misc : Blank
 ALS Vial : 2 Sample Multiplier: 1

Quant Time: Aug 06 11:32:06 2014
 Quant Method : C:\msdchem\1\METHODS (b) (7) RP1581_051914.M
 Quant Title : Freedom Industries (C)
 QLast Update : Wed Aug 06 11:28:36 2014
 Response via : Initial Calibration



Quantitation Report (QT Reviewed)

Data Path : C:\msdchem\1\DATA\ (b) (7) \RP1581\
Data File : 051914_003.D
Acq On : 19 May 2014 12:55 pm
Operator : (b) (7)(C)
Sample : 0.5 ug/mL
Misc : 0.5
ALS Vial : 3 Sample Multiplier: 1

Quant Time: Aug 06 11:40:07 2014
Quant Method : C:\msdchem\1\METHODS\ (b) (7) \RP1581_051914.M
Quant Title : Freedom Industries
QLast Update : Wed Aug 06 11:32:39 2014
Response via : Initial Calibration

Compound	(b) (7)(C)	Response	Conc	Units	Dev (Min)


System Monitoring Compounds					
1) Cyclohexanol	4.613	57	156539	48.53 ug/mL	0.00
Target Compounds					
2) MCHM (2 peaks combined)	6.530	55	842m	2.86 ug/mL	Qvalue
3) PPH	7.784	94	548m	4.27 ug/mL	
4) 1,4-CHDM	8.592	95	280m	4.08 ug/mL	
5) DMCH-1,4-DC	9.152	81	646m	3.18 ug/mL	
6) di-PPH	9.732	59	401m	2.77 ug/mL	

(#) = qualifier out of range (m) = manual integration (+) = signals summed

Quantitation Report (QT Reviewed)

Data Path : C:\msdchem\1\DATA (b) (7)(C) RP1581\
 Data File : 051914_003.D
 Acq On : 19 May 2014 12:55 pm
 Operator : (b) (7)(C)
 Sample : 0.5 ug/mL
 Misc : 0.5
 ALS Vial : 3 Sample Multiplier: 1

Quant Time: Aug 06 11:40:07 2014
 Quant Method : C:\msdchem\1\METHODS (b) (7) RP1581_051914.M
 Quant Title : Freedom Industries
 QLast Update : Wed Aug 06 11:32:39 2014
 Response via : Initial Calibration



Quantitation Report (QT Reviewed)

Data Path : C:\msdchem\1\DATA\ (b) (7) \RP1581\
Data File : 051914_004.D
Acq On : 19 May 2014 1:28 pm
Operator : (b) (7)
Sample : 1 ug/mL
Misc : 1
ALS Vial : 4 Sample Multiplier: 1

Quant Time: Aug 06 11:41:21 2014
Quant Method : C:\msdchem\1\METHODS (b) (7) RP1581_051914.M
Quant Title : Freedom Industries (b) (7)
QLast Update : Wed Aug 06 11:32:39 2014
Response via : Initial Calibration

Compound	(b) (7)(C)		Response	Conc	Units	Dev (Min)


System Monitoring Compounds						
1) Cyclohexanol	4.613	57	155813	48.30	ug/mL	0.00
Target Compounds						
2) MCHM (2 peaks combined)	6.530	55	2156m	3.18	ug/mL	Qvalue
3) PPH	7.784	94	1506m	4.45	ug/mL	
4) 1,4-CHDM	8.582	95	545m	4.20	ug/mL	
5) DMCH-1,4-DC	9.152	81	1398m	3.41	ug/mL	
6) di-PPH	9.846	59	1017m	3.04	ug/mL	

(#) = qualifier out of range (m) = manual integration (+) = signals summed

Quantitation Report (QT Reviewed)

Data Path : C:\msdchem\1\DATA (b) (7)(C) RP1581\
 Data File : 051914_004.D
 Acq On : 19 May 2014 1:28 pm
 Operator : (b) (7)(C)
 Sample : 1 ug/mL
 Misc : 1
 ALS Vial : 4 Sample Multiplier: 1

Quant Time: Aug 06 11:41:21 2014
 Quant Method : C:\msdchem\1\METHODS (b) (7) RP1581_051914.M
 Quant Title : Freedom Industries (C)
 QLast Update : Wed Aug 06 11:32:39 2014
 Response via : Initial Calibration



Quantitation Report (QT Reviewed)

Data Path : C:\msdchem\1\DATA (b) (7)(C) RP1581\
Data File : 051914_005.D
Acq On : 19 May 2014 2:00 pm
Operator : (b) (7)
Sample : (C) ug/mL
Misc : 3
ALS Vial : 5 Sample Multiplier: 1

Quant Time: Aug 06 11:42:27 2014
Quant Method : C:\msdchem\1\METHODS (b) (7)(C) RP1581_051914.M
Quant Title : Freedom Industries
QLast Update : Wed Aug 06 11:32:39 2014
Response via : Initial Calibration

Compound	(b) (7)(C)		Response	Conc	Units	Dev (Min)


System Monitoring Compounds						
1) Cyclohexanol	4.623	57	169255	52.47	ug/mL	0.01
Target Compounds						
2) MCHM (2 peaks combined)	6.530	55	7148m	4.37	ug/mL	Qvalue
3) PPH	7.774	94	6493m	5.38	ug/mL	
4) 1,4-CHDM	8.582	95	2435m	5.07	ug/mL	
5) DMCH-1,4-DC	9.152	81	5114m	4.51	ug/mL	
6) di-PPH	9.846	59	3031m	3.94	ug/mL	

(#) = qualifier out of range (m) = manual integration (+) = signals summed

Quantitation Report (QT Reviewed)

Data Path : C:\msdchem\1\DATA\ (b) (7)(C) RP1581\
 Data File : 051914_005.D
 Acq On : 19 May 2014 2:00 pm
 Operator : (b) (7)
 Sample : 3 ug/mL
 Misc : 3
 ALS Vial : 5 Sample Multiplier: 1

Quant Time: Aug 06 11:42:27 2014
 Quant Method : C:\msdchem\1\METHODS\ (b) (7)(C) RP1581_051914.M
 Quant Title : Freedom Industries
 QLast Update : Wed Aug 06 11:32:39 2014
 Response via : Initial Calibration



Data Path : C:\msdchem\1\DATA(b) (7)(C)\RP1581\
 Data File : 051914_006.D
 Acq On : 19 May 2014 2:32 pm
 Operator : (b) (7)
 Sample : 5 ug/ml
 Misc : 5
 ALS Vial : 6 Sample Multiplier: 1

Quant Time: Aug 06 11:43:33 2014
 Quant Method : C:\msdchem\1\METHODS(b) (7)(C)\RP1581_051914.M
 Quant Title : Freedom Industries
 QLast Update : Wed Aug 06 11:32:39 2014
 Response via : Initial Calibration

Compound	R.T.	QIon	Response	Conc	Units	Dev (Min)


System Monitoring Compounds						
1) Cyclohexanol	4.623	57	157595	48.86	ug/mL	0.01
Target Compounds						
2) MCHM (2 peaks combined)	6.530	55	10861m	5.26	ug/mL	Qvalue
3) PPH	7.774	94	10474m	6.12	ug/mL	
4) 1,4-CHDM	8.582	95	4360m	5.95	ug/mL	
5) DMCH-1,4-DC	9.152	81	7673m	5.27	ug/mL	
6) di-PPH	9.846	59	4799m	4.73	ug/mL	

(#) = qualifier out of range (m) = manual integration (+) = signals summed

Quantitation Report (QT Reviewed)

Data Path : C:\msdchem\1\DATA\ (b) (7)(C) RP1581\
 Data File : 051914_006.D
 Acq On : 19 May 2014 2:32 pm
 Operator : (b) (7)(C)
 Sample : 5 ug/ml
 Misc : 5
 ALS Vial : 6 Sample Multiplier: 1

Quant Time: Aug 06 11:43:33 2014
 Quant Method : C:\msdchem\1\METHODS\ (b) (7)(C) P1581_051914.M
 Quant Title : Freedom Industries
 QLast Update : Wed Aug 06 11:32:39 2014
 Response via : Initial Calibration



Data Path : C:\msdchem\1\DATA (b) (7)(C) \RP1581\
Data File : 051914_007.D
Acq On : 19 May 2014 3:04 pm
Operator : (b) (7)
Sample : 10 ug/mL
Misc : 10
ALS Vial : 7 Sample Multiplier: 1

Quant Time: Aug 06 11:44:19 2014
Quant Method : C:\msdchem\1\METHODS\ (b) (7) \RP1581_051914.M
Quant Title : Freedom Industries
QLast Update : Wed Aug 06 11:32:39 2014
Response via : Initial Calibration

Compound	R.T.	QIon	Response	Conc	Units	Dev (Min)


System Monitoring Compounds						
1) Cyclohexanol	4.623	57	151352	46.92	ug/mL	0.01
Target Compounds						Qvalue
2) MCHM (2 peaks combined)	6.530	55	22258m	8.00	ug/mL	
3) PPH	7.774	94	21573	8.19	ug/mL#	55
4) 1,4-CHDM	8.582	95	8585m	7.90	ug/mL	
5) DMCH-1,4-DC	9.152	81	18423m	8.46	ug/mL	
6) di-PPH	9.846	59	11129m	7.55	ug/mL	

(#) = qualifier out of range (m) = manual integration (+) = signals summed

Quantitation Report (QT Reviewed)

Data Path : C:\msdchem\1\DATA (b) (7)(C) RP1581\
 Data File : 051914_007.D
 Acq On : 19 May 2014 3:04 pm
 Operator : (b) (7)(C)
 Sample : 10 ug/mL
 Misc : 10
 ALS Vial : 7 Sample Multiplier: 1

Quant Time: Aug 06 11:44:19 2014
 Quant Method : C:\msdchem\1\METHODS (b) (7)(C) RP1581_051914.M
 Quant Title : Freedom Industries
 QLast Update : Wed Aug 06 11:32:39 2014
 Response via : Initial Calibration



Data Path : C:\msdchem\1\DATA\ (b) (7)(C) RP1581\
Data File : 051914_008.D
Acq On : 19 May 2014 3:36 pm
Operator : (b) (7)
Sample : 50 ug/mL
Misc : 50
ALS Vial : 8 Sample Multiplier: 1

Quant Time: Aug 06 11:45:25 2014
Quant Method : C:\msdchem\1\METHODS\AMH_RP1581_051914.M
Quant Title : Freedom Industries
QLast Update : Wed Aug 06 11:32:39 2014
Response via : Initial Calibration

Compound	(b) (7)(C)	Response	Conc	Units	Dev (Min)


System Monitoring Compounds					
1) Cyclohexanol	4.613	57	168333	52.19 ug/mL	0.00
Target Compounds					
2) MCHM (2 peaks combined)	6.530	55	163896m	41.96 ug/mL	Qvalue
3) PPH	7.773	94	189499	39.50 ug/mL	94
4) 1,4-CHDM	8.582	95	74607m	38.26 ug/mL	
5) DMCH-1,4-DC	9.152	81	130770	41.79 ug/mL#	91
6) di-PPH	9.846	59	83192m	39.69 ug/mL	

(#) = qualifier out of range (m) = manual integration (+) = signals summed

Quantitation Report (QT Reviewed)

Data Path : C:\msdchem\1\DATA (b) (7)(C) RP1581\
 Data File : 051914_008.D
 Acq On : 19 May 2014 3:36 pm
 Operator : (b) (7)(C)
 Sample : 50 ug/mL
 Misc : 50
 ALS Vial : 8 Sample Multiplier: 1

Quant Time: Aug 06 11:45:25 2014
 Quant Method : C:\msdchem\1\METHODS (b) (7)(C) RP1581_051914.M
 Quant Title : Freedom Industries
 QLast Update : Wed Aug 06 11:32:39 2014
 Response via : Initial Calibration



Data Path : C:\msdchem\1\DATA\ (b) (7)(C) RP1581\
Data File : 051914_009.D
Acq On : 19 May 2014 4:09 pm
Operator : (b) (7)(C)
Sample : 100 ug/mL
Misc : 100
ALS Vial : 9 Sample Multiplier: 1

Quant Time: Aug 06 11:46:47 2014
Quant Method : C:\msdchem\1\METHODS\ (b) (7)(C) P1581_051914.M
Quant Title : Freedom Industries
QLast Update : Wed Aug 06 11:32:39 2014
Response via : Initial Calibration

Compound	R.T.	QIon	Response	Conc	Units	Dev (Min)


System Monitoring Compounds						
1) Cyclohexanol	4.623	57	167242	51.85	ug/mL	0.01
Target Compounds						
2) MCHM (2 peaks combined)	6.530	55	367546m	90.81	ug/mL	Qvalue
3) PPH	7.774	94	441478	86.47	ug/mL	95
4) 1,4-CHDM	8.582	95	186886m	89.89	ug/mL	
5) DMCH-1,4-DC	9.152	81	293609	90.11	ug/mL#	91
6) di-PPH	9.846	59	214165m	98.11	ug/mL	

(#) = qualifier out of range (m) = manual integration (+) = signals summed

Quantitation Report (QT Reviewed)

Data Path : C:\msdchem\1\DATA (b) (7)(C) RP1581\
 Data File : 051914_009.D
 Acq On : 19 May 2014 4:09 pm
 Operator : (b) (7)(C)
 Sample : 100 ug/mL
 Misc : 100
 ALS Vial : 9 Sample Multiplier: 1

Quant Time: Aug 06 11:46:47 2014
 Quant Method : C:\msdchem\1\METHODS (b) (7)(C) RP1581_051914.M
 Quant Title : Freedom Industries
 QLast Update : Wed Aug 06 11:32:39 2014
 Response via : Initial Calibration



Data Path : C:\msdchem\1\DATA (b) (7)(C) \RP1581\
Data File : 051914_010.D
Acq On : 19 May 2014 4:42 pm
Operator : (b) (7)(C)
Sample : g/mL
Misc : 200
ALS Vial : 10 Sample Multiplier: 1

Quant Time: Aug 06 11:47:37 2014 (b) (7)
Quant Method : C:\msdchem\1\METHODS (C) RP1581_051914.M
Quant Title : Freedom Industries
QLast Update : Wed Aug 06 11:32:39 2014
Response via : Initial Calibration

Compound	(b) (7)(C)	Response	Conc	Units	Dev (Min)


System Monitoring Compounds					
1) Cyclohexanol	4.623	57	163415	50.66 ug/mL	0.01
Target Compounds					
2) MCHM (2 peaks combined)	6.530	55	789712m	192.05 ug/mL	Qvalue
3) PPH	7.784	94	1023617	194.99 ug/mL	96
4) 1,4-CHDM	8.582	95	410254m	192.61 ug/mL	
5) DMCH-1,4-DC	9.162	81	639126	192.62 ug/mL#	92
6) di-PPH	9.846	59	527324m	237.79 ug/mL	

(#) = qualifier out of range (m) = manual integration (+) = signals summed

Quantitation Report (QT Reviewed)

Data Path : C:\msdchem\1\DATA (b) (7)(C) RP1581\
 Data File : 051914_010.D
 Acq On : 19 May 2014 4:42 pm
 Operator : (b) (7)(C)
 Sample : 200 ug/mL
 Misc : 200
 ALS Vial : 10 Sample Multiplier: 1

Quant Time: Aug 06 11:47:37 2014
 Quant Method : C:\msdchem\1\METHODS (b) (7)(C) RP1581_051914.M
 Quant Title : Freedom Industries
 QLast Update : Wed Aug 06 11:32:39 2014
 Response via : Initial Calibration



METHOD CONTROL PARAMETERS

Method Information For: C:\MSDCHEM\1\METHODS\ (b) (7) (C) RP1581.M
Method Sections To Run:

() Save Copy of Method With Data
() Instrument Control Pre-Run Cmd/Macro =
() Data Analysis Pre-Run Cmd/Macro =
(X) Data Acquisition
(X) Data Analysis
() Instrument Control Post-Run Cmd/Macro =
() Data Analysis Post-Run Cmd/Macro =
Method Comments:

END OF METHOD CONTROL PARAMETERS

INSTRUMENT CONTROL PARAMETERS: Instrument 1

C:\MSDCHEM\1\METHODS(b) (7) RP1581.M
Thu Jun 05 14:28:38 2014

Control Information

Sample Inlet : GC
Injection Source : GC ALS
Mass Spectrometer : Enabled

Oven
Equilibration Time 0.5 min
Oven Program On
35 °C for 2 min
then 20 °C/min to 325 °C for 3 min
Run Time 19.5 min

Front Injector
Syringe Size 10 µL
Injection Volume 1 µL
Injection Repetitions 1
Solvent A Washes (PreInj) 2
Solvent A Washes (PostInj) 3
Solvent A Volume 8 µL
Solvent B Washes (PreInj) 2
Solvent B Washes (PostInj) 3
Solvent B Volume 8 µL
Sample Washes 1
Sample Wash Volume 8 µL
Sample Pumps 3
Dwell Time (PreInj) 0 min
Dwell Time (PostInj) 0 min
Solvent Wash Draw Speed 300 µL/min
Solvent Wash Dispense Speed 6000 µL/min
Sample Wash Draw Speed 300 µL/min
Sample Wash Dispense Speed 6000 µL/min
Injection Dispense Speed 6000 µL/min
Viscosity Delay 0 sec
Sample Depth Disabled

Back Injector

Front SS Inlet He
Mode Split
Heater On 275 °C
Pressure On 16.494 psi
Total Flow On 47.305 mL/min
Septum Purge Flow On 3 mL/min
Gas Saver On 30 mL/min After 2 min
Split Ratio 20 :1
Split Flow 42.195 mL/min

Back SS Inlet He
Mode Split
Heater Off
Pressure Off
Total Flow Off
Septum Purge Flow Off
Gas Saver Off
Split Ratio 100 :1

Split Flow

0 mL/min

Thermal Aux 2 {MSD Transfer Line}

Heater On

Temperature Program On

280 °C for 0 min

Run Time 19.5 min

Column #1

DB-5ms: 1653.60651

DB-5ms

325 °C: 30 m x 250 µm x 0.25 µm

In: Front SS Inlet He

Out: Vacuum

(Initial)

35 °C

Pressure 16.494 psi

Flow 2.1098 mL/min

Average Velocity 52.535 cm/sec

Holdup Time 0.95175 min

Flow Program Off

2.1098 mL/min for 0 min

Run Time 19.5 min

Front Detector µECD

Heater On 300 °C

Anode Flow Off

Makeup Flow On 20 mL/min

Const Col + Makeup Off

Electrometer Off

Aux Pressure 1 He

Pressure Program Off

10 psi for 0 min

Run Time 19.5 min

Aux Pressure 2 He

Pressure Program Off

10 psi for 0 min

Run Time 19.5 min

Aux Pressure 3 He

Pressure Program Off

10 psi for 0 min

Run Time 19.5 min

Signals

Front Signal Save Off

Back Inlet Save Off

Test Plot Save Off

Test Plot Save Off

MS ACQUISITION PARAMETERS

* General Information

Tune File : dftppnew051614.u
Acquisition Mode : Scan

(b) Information
(7) -----

Solvent Delay : 3.50 min
EMV Mode : Relative
Relative Voltage : -71
Resulting EM Voltage : 1812

[Scan Parameters]

Low Mass : 35.0
High Mass : 500.0
Threshold : 50
Sample # : 3 A/D Samples 8
Plot 2 low mass : 50.0
Plot 2 high mass : 550.0

[MSZones]

MS Source : 230 C maximum 250 C
MS Quad : 150 C maximum 200 C

END OF MS ACQUISITION PARAMETERS

TUNE PARAMETERS for SN: US73347311

Trace Ion Detection is OFF.

EMISSION : 34.610
ENERGY : 69.922
REPELLER : 22.249
IONFOCUS : 90.157
ENTRANCE_LE : 0.000
EMVOLTS : 1882.353

Actual EMV : 1811.77
GAIN FACTOR : 2.25

AMUGAIN : 2075.000
AMUOFFSET : 126.750
FILAMENT : 1.000
DCPOLARITY : 0.000
ENTLENSOFFS : 17.318@ 3 17.318@ 50 15.310@ 69 13.302@131 13.804
@219 14.055@414 14.055@502 14.055@1049
MASSGAIN : -453.000
MASSOFFSET : -39.000

END OF TUNE PARAMETERS

END OF INSTRUMENT CONTROL PARAMETERS

DATA ANALYSIS PARAMETERS

Method Name: C:\MSDCHEM\1\METHODS\RP1581.M

Percent Report Settings

Sort By: Signal

Output Destination

Screen: No
Printer: No
File: No

Integration Events: Meth Default

Generate Report During Run Method: Yes

Signal Correlation Window: 0.020

Qualitative Report Settings

Peak Location of Unknown: Apex

Library to Search Minimum Quality
C:\Database\NIST05a.L 0

Integration Events: Meth Default

Report Type: Summary

Output Destination

Screen: No
Printer: No
File: No

Generate Report During Run Method: Yes

Quantitative Report Settings

Report Type: Detailed

Output Destination

Screen: No
Printer: No
File: epatemp.txt

Generate Report During Run Method: No

Freedom Industries

Calibration Last Updated: Tue Jun 03 14:32:34 2014

Reference Window: 2.00 Minutes

Non-Reference Window: 1.00 Minutes

AMH_RP1581.M Thu Jun 05 14:28:37 2014

Correlation Window: 0.10 minutes
Default Multiplier: 1.00
Default Sample Concentration: 0.00

Compound Information

1) Cyclohexanol

()

Ret. Time 4.613 min., Extract & Integrate from 4.113 to 5.113 min.

Signal	Rel Resp.	Pct. Unc.(rel)	Integration
Tgt 57.10			*** METH DEFAULT ***
Q1 82.10	57.10	20.0	*** METH DEFAULT ***
Q2 67.10	43.80	20.0	*** METH DEFAULT ***

Lvl ID	Conc (ug/mL)	Response
1	50.000	156539
2	50.000	155813
3	50.000	169255
4	50.000	157595
5	50.000	151352
6	50.000	168333
7	50.000	167242
8	50.000	163415

Qualifier Peak Analysis ON
Curve Fit: Avg. RF

2) MCHM (2 peaks combined)

()

Ret. Time 6.530 min., Extract & Integrate from 6.030 to 7.030 min.

Signal	Rel Resp.	Pct. Unc.(rel)	Integration
Tgt 55.10			*** METH DEFAULT ***
Q1 97.10	39.40	20.0	*** METH DEFAULT ***
Q2 81.10	30.30	20.0	*** METH DEFAULT ***

Lvl ID	Conc (ug/mL)	Response
1	0.500	653
2	1.000	1464
3	3.000	4958
4	5.000	7581
5	10.000	15368
6	50.000	118295
7	100.000	265510
8	200.000	550216

Qualifier Peak Analysis ON
Curve Fit: Linear

3) PPH

()

Ret. Time 7.784 min., Extract & Integrate from 7.284 to 8.284 min.

Signal	Rel Resp.	Pct. Unc.(rel)	Integration
Tgt 94.10			*** METH DEFAULT ***
Q1 77.10	23.70	20.0	*** METH DEFAULT ***
Q2 152.10	18.60	20.0	*** METH DEFAULT ***

Lvl ID	Conc (ug/mL)	Response
1	0.500	548
2	1.000	1506
3	3.000	6257
4	5.000	10053
5	10.000	21573
6	50.000	189499
7	100.000	441478
8	200.000	1023617

Qualifier Peak Analysis ON
Curve Fit: Linear

4) 1,4-CHDM

()

Ret. Time 8.592 min., Extract & Integrate from 8.092 to 9.092 min.

Signal	Rel Resp.	Pct. Unc.(rel)	Integration
Tgt 95.10			*** METH DEFAULT ***
Q1 67.10	41.40	20.0	*** METH DEFAULT ***
Q2 41.10	53.20	20.0	*** METH DEFAULT ***

Lvl ID	Conc (ug/mL)	Response
1	0.500	280
2	1.000	545
3	3.000	1859
4	5.000	3449
5	10.000	6555
6	50.000	58043
7	100.000	147139
8	200.000	315052

Qualifier Peak Analysis ON
Curve Fit: Linear

5) DMCH-1,4-DC

()

Ret. Time 9.152 min., Extract & Integrate from 8.652 to 9.652 min.

Signal	Rel Resp.	Pct. Unc.(rel)	Integration
Tgt 81.10			*** METH DEFAULT ***
Q1 108.10	36.20	20.0	*** METH DEFAULT ***
Q2 140.10	63.60	20.0	*** METH DEFAULT ***

Lvl ID	Conc (ug/mL)	Response
1	0.500	646
2	1.000	1399
3	3.000	5113
4	5.000	7673
5	10.000	17864
6	50.000	130770
7	100.000	293609
8	200.000	639126

Qualifier Peak Analysis ON
Curve Fit: Linear

6) di-PPH

()

Ret. Time 9.856 min., Extract & Integrate from 9.356 to 10.356 min.

Signal	Rel Resp.	Pct. Unc.(rel)	Integration
Tgt 59.10			*** METH DEFAULT ***
Q1 94.10	30.00	20.0	*** METH DEFAULT ***
Q2 210.20	9.00	20.0	*** METH DEFAULT ***

Lvl ID	Conc (ug/mL)	Response
1	0.500	221
2	1.000	726
3	3.000	2147
4	5.000	3145
5	10.000	7351
6	50.000	58083
7	100.000	154086
8	200.000	375689

Qualifier Peak Analysis ON
Curve Fit: Linear

END OF DATA ANALYSIS PARAMETERS

Thu Jun 05 14:28:39 2014