



# Analytical Laboratory

Analytical Lab  
Page 1 of 28

13339 Hagers Ferry Road  
Huntersville, NC 28078-7929  
McGuire Nuclear Complex - MG03A2  
Phone: 980-875-5245 Fax: 980-875-4349

## Order Summary Report

**Order Number:** J11120115

**Customer Name(s):** Bill Kennedy, Melonie Martin, Wayne Chapman, Tom Johnson

**Customer Address:** 3195 Pine Hall Rd  
Mailcode: Belews Steam Station  
Belews Creek, NC 28012

**Lab Contact:** Jason C Perkins **Phone:** 980-875-5348

**Report Authorized By:** \_\_\_\_\_ **Date:** 1/4/2012  
**(Signature)**

---

### Program Comments:

Please contact the Program Manager (Jason C Perkins) with any questions regarding this report.

### Data Flags & Calculations:

Any analytical tests or individual analytes within a test flagged with a Qualifier indicate a deviation from the method quality system or quality control requirement. The qualifier description is found at the end of the Certificate of Analysis (sample results) under the qualifiers heading. All results are reported on a dry weight basis unless otherwise noted.

### Data Package:

This data package includes analytical results that are applicable only to the samples described in this narrative. An estimation of the uncertainty of measurement for the results in the report is available upon request. This report shall not be reproduced, except in full, without the written consent of the Analytical Laboratory. Please contact the Analytical laboratory with any questions. The order of individual sections within this report is as follows:

*Job Summary Report, Sample Identification, Technical Validation of Data Package, Analytical Laboratory Certificate of Analysis, Analytical Laboratory QC Reports, Sub-contracted Laboratory Results, Customer Specific Data Sheets, Reports & Documentation, Customer Database Entries, Test Case Narratives, Chain of Custody (COC)*

### Certification:

The Analytical Laboratory holds the following State Certifications : North Carolina (DENR) Certificate #248, South Carolina (DHEC) Laboratory ID # 99005. Contact the Analytical Laboratory for definitive information about the certification status of specific methods.

## Sample ID's & Descriptions:

Sample ID	Plant/Station	Collection Date and Time	Collected By	Sample Description
2011026219	BELEWS	14-Dec-11 11:30 AM	DAVID MORRIS	FGD Purge Eff
2011026342	BELEWS	14-Dec-11 11:30 AM	DAVID MORRIS	EQ TANK EFF.
2011026343	BELEWS	14-Dec-11 11:30 AM	DAVID MORRIS	BIOREACTOR 1 INF.
2011026344	BELEWS	14-Dec-11 11:30 AM	DAVID MORRIS	BIOREACTOR 2 INF.
2011026345	BELEWS	14-Dec-11 11:30 AM	DAVID MORRIS	BIOREACTOR 2 EFF.
2011026346	BELEWS	14-Dec-11 11:30 AM	DAVID MORRIS	FILTER BLANK
2011026347	BELEWS	14-Dec-11 11:30 AM	DAVID MORRIS	Trip Blank
2011026356	BELEWS	14-Dec-11 11:40 AM	DAVID MORRIS	BIOREACTOR 1 INF.
2011026357	BELEWS	14-Dec-11 11:40 AM	DAVID MORRIS	HG BLANK BIOREACTOR 1 INF.
2011026358	BELEWS	14-Dec-11 11:50 AM	DAVID MORRIS	BIOREACTOR 2 INF.
2011026359	BELEWS	14-Dec-11 11:50 AM	DAVID MORRIS	Hg Blk BioReactor 2 Inf
2011026360	BELEWS	14-Dec-11 11:45 AM	DAVID MORRIS	BIOREACTOR 2 EFF.
2011026361	BELEWS	14-Dec-11 11:45 AM	DAVID MORRIS	Hg Blk BioReactor 2 Eff
13 Total Samples				

# Technical Validation Review

## Checklist:

COC and .pdf report are in agreement with sample totals and analyses (compliance programs and procedures).

☒ Yes

☐ No

All Results are less than the laboratory reporting limits.

☐ Yes

☒ No

All laboratory QA/QC requirements are acceptable.

☒ Yes

☐ No

The Vendor Laboratories have been qualified by the Analytical Laboratory

Yes

## Report Sections Included:

☒ Job Summary Report

☒ Sample Identification

☒ Technical Validation of Data Package

☒ Analytical Laboratory Certificate of Analysis

☐ Analytical Laboratory QC Report

☒ Sub-contracted Laboratory Results

☐ Customer Specific Data Sheets, Reports, & Documentation

☐ Customer Database Entries

☒ Chain of Custody

☒ Electronic Data Deliverable (EDD) Sent Separately

Reviewed By: Mary Ann Ogle

Date: 1/4/2012

**Certificate of Laboratory Analysis***This report shall not be reproduced, except in full.***Order # J11120115**Site: FGD Purge Eff  
Collection Date: 14-Dec-11 11:30 AMSample #: **2011026219**  
Matrix: OTHER

Analyte	Result	Units	Qualifiers	RDL	Method	Analysis Date/Time	Analyst
<b><u>INORGANIC IONS BY IC</u></b>							
Bromide	84	mg/L		5	EPA 300.0	19-Dec-11 19:30	JAHERMA
<b><u>MERCURY (COLD VAPOR) IN WATER</u></b>							
Mercury (Hg)	229	ug/L		5	EPA 245.1	16-Dec-11 09:09	AGIBBS
<b><u>TOTAL RECOVERABLE METALS BY ICP</u></b>							
Boron (B)	146	mg/L		0.5	EPA 200.7	16-Dec-11 10:17	DJSULL1
<b><u>DISSOLVED METALS BY ICP-MS</u></b>							
Selenium (Se)	419	ug/L		10	EPA 200.8	20-Dec-11 10:55	MHH7131
<b><u>TOTAL RECOVERABLE METALS BY ICP-MS</u></b>							
Arsenic (As)	170	ug/L		10	EPA 200.8	21-Dec-11 09:43	MHH7131
Chromium (Cr)	172	ug/L		10	EPA 200.8	21-Dec-11 09:43	MHH7131
Copper (Cu)	108	ug/L		10	EPA 200.8	21-Dec-11 09:43	MHH7131
Nickel (Ni)	166	ug/L		10	EPA 200.8	21-Dec-11 09:43	MHH7131
Selenium (Se)	4440	ug/L		10	EPA 200.8	21-Dec-11 09:43	MHH7131
Silver (Ag)	< 10.0	ug/L		10	EPA 200.8	21-Dec-11 09:43	MHH7131
Zinc (Zn)	181	ug/L		20	EPA 200.8	21-Dec-11 09:43	MHH7131
<b><u>SELENIUM SPECIATION</u></b>							
Vendor Parameter	Complete				V_AS&C		
<b><u>TOTAL DISSOLVED SOLIDS</u></b>							
TDS	15000	mg/L		200	SM2540C	19-Dec-11 14:36	AGIBBS

Site: EQ TANK EFF.  
Collection Date: 14-Dec-11 11:30 AMSample #: **2011026342**  
Matrix: OTHER

Analyte	Result	Units	Qualifiers	RDL	Method	Analysis Date/Time	Analyst
<b><u>MERCURY (COLD VAPOR) IN WATER</u></b>							
Mercury (Hg)	208	ug/L		2.5	EPA 245.1	16-Dec-11 09:12	AGIBBS
<b><u>TOTAL RECOVERABLE METALS BY ICP</u></b>							
Boron (B)	152	mg/L		0.5	EPA 200.7	16-Dec-11 10:21	DJSULL1
<b><u>DISSOLVED METALS BY ICP-MS</u></b>							
Selenium (Se)	180	ug/L		10	EPA 200.8	20-Dec-11 10:55	MHH7131

# Certificate of Laboratory Analysis

Analytical Lab  
Page 5 of 28

*This report shall not be reproduced, except in full.*

Order # J11120115

Site: EQ TANK EFF.

Collection Date: 14-Dec-11 11:30 AM

Sample #: 2011026342

Matrix: OTHER

Analyte	Result	Units	Qualifiers	RDL	Method	Analysis Date/Time	Analyst
<b>TOTAL RECOVERABLE METALS BY ICP-MS</b>							
Arsenic (As)	164	ug/L		10	EPA 200.8	21-Dec-11 09:46	MHH7131
Chromium (Cr)	169	ug/L		10	EPA 200.8	21-Dec-11 09:46	MHH7131
Copper (Cu)	106	ug/L		10	EPA 200.8	21-Dec-11 09:46	MHH7131
Nickel (Ni)	162	ug/L		10	EPA 200.8	21-Dec-11 09:46	MHH7131
Selenium (Se)	4160	ug/L		10	EPA 200.8	21-Dec-11 09:46	MHH7131
Silver (Ag)	< 10.0	ug/L		10	EPA 200.8	21-Dec-11 09:46	MHH7131
Zinc (Zn)	178	ug/L		20	EPA 200.8	21-Dec-11 09:46	MHH7131

Site: BIOREACTOR 1 INF.

Collection Date: 14-Dec-11 11:30 AM

Sample #: 2011026343

Matrix: OTHER

Analyte	Result	Units	Qualifiers	RDL	Method	Analysis Date/Time	Analyst
<b>TOTAL RECOVERABLE METALS BY ICP</b>							
Boron (B)	146	mg/L		0.5	EPA 200.7	16-Dec-11 10:25	DJSULL1
<b>DISSOLVED METALS BY ICP-MS</b>							
Selenium (Se)	141	ug/L		10	EPA 200.8	20-Dec-11 10:55	MHH7131
<b>TOTAL RECOVERABLE METALS BY ICP-MS</b>							
Arsenic (As)	< 10.0	ug/L		10	EPA 200.8	21-Dec-11 09:49	MHH7131
Chromium (Cr)	< 10.0	ug/L		10	EPA 200.8	21-Dec-11 09:49	MHH7131
Copper (Cu)	< 10.0	ug/L		10	EPA 200.8	21-Dec-11 09:49	MHH7131
Nickel (Ni)	39.7	ug/L		10	EPA 200.8	21-Dec-11 09:49	MHH7131
Selenium (Se)	133	ug/L		10	EPA 200.8	21-Dec-11 09:49	MHH7131
Silver (Ag)	< 10.0	ug/L		10	EPA 200.8	21-Dec-11 09:49	MHH7131
Zinc (Zn)	< 20.0	ug/L		20	EPA 200.8	21-Dec-11 09:49	MHH7131

## SELENIUM SPECIATION

Vendor Parameter Complete V\_AS&C

Site: BIOREACTOR 2 INF.

Collection Date: 14-Dec-11 11:30 AM

Sample #: 2011026344

Matrix: OTHER

Analyte	Result	Units	Qualifiers	RDL	Method	Analysis Date/Time	Analyst
<b>TOTAL RECOVERABLE METALS BY ICP</b>							
Boron (B)	138	mg/L		0.5	EPA 200.7	16-Dec-11 10:29	DJSULL1

# Certificate of Laboratory Analysis

Analytical Lab  
Page 6 of 28

*This report shall not be reproduced, except in full.*

Order # J11120115

Site: BIOREACTOR 2 INF.

Collection Date: 14-Dec-11 11:30 AM

Sample #: 2011026344

Matrix: OTHER

Analyte	Result	Units	Qualifiers	RDL	Method	Analysis Date/Time	Analyst
<b>TOTAL RECOVERABLE METALS BY ICP-MS</b>							
Arsenic (As)	< 10.0	ug/L		10	EPA 200.8	21-Dec-11 09:52	MHH7131
Chromium (Cr)	< 10.0	ug/L		10	EPA 200.8	21-Dec-11 09:52	MHH7131
Copper (Cu)	< 10.0	ug/L		10	EPA 200.8	21-Dec-11 09:52	MHH7131
Nickel (Ni)	< 10.0	ug/L		10	EPA 200.8	21-Dec-11 09:52	MHH7131
Selenium (Se)	14.5	ug/L		10	EPA 200.8	21-Dec-11 09:52	MHH7131
Silver (Ag)	< 10.0	ug/L		10	EPA 200.8	21-Dec-11 09:52	MHH7131
Zinc (Zn)	< 20.0	ug/L		20	EPA 200.8	21-Dec-11 09:52	MHH7131

Site: BIOREACTOR 2 EFF.

Collection Date: 14-Dec-11 11:30 AM

Sample #: 2011026345

Matrix: OTHER

Analyte	Result	Units	Qualifiers	RDL	Method	Analysis Date/Time	Analyst
<b>INORGANIC IONS BY IC</b>							
Bromide	82	mg/L		5	EPA 300.0	19-Dec-11 19:46	JAHERMA
<b>MERCURY (COLD VAPOR) IN WATER</b>							
Mercury (Hg)	< 1.00	ug/L		1	EPA 245.1	16-Dec-11 09:14	AGIBBS
<b>TOTAL RECOVERABLE METALS BY ICP</b>							
Boron (B)	129	mg/L		0.5	EPA 200.7	16-Dec-11 10:32	DJSULL1
<b>TOTAL RECOVERABLE METALS BY ICP-MS</b>							
Arsenic (As)	< 5.00	ug/L		5	EPA 200.8	21-Dec-11 09:55	MHH7131
Chromium (Cr)	< 5.00	ug/L		5	EPA 200.8	21-Dec-11 09:55	MHH7131
Copper (Cu)	< 5.00	ug/L		5	EPA 200.8	21-Dec-11 09:55	MHH7131
Nickel (Ni)	< 5.00	ug/L		5	EPA 200.8	21-Dec-11 09:55	MHH7131
Selenium (Se)	< 5.00	ug/L		5	EPA 200.8	21-Dec-11 09:55	MHH7131
Silver (Ag)	< 5.00	ug/L		5	EPA 200.8	21-Dec-11 09:55	MHH7131
Zinc (Zn)	< 10.0	ug/L		10	EPA 200.8	21-Dec-11 09:55	MHH7131
<b>SELENIUM SPECIATION</b>							
Vendor Parameter	Complete				V_AS&C		

Site: FILTER BLANK

Collection Date: 14-Dec-11 11:30 AM

Sample #: 2011026346

Matrix: OTHER

Analyte	Result	Units	Qualifiers	RDL	Method	Analysis Date/Time	Analyst
<b>DISSOLVED METALS BY ICP-MS</b>							
Selenium (Se)	< 1.000	ug/L		1	EPA 200.8	20-Dec-11 10:55	MHH7131

# Certificate of Laboratory Analysis

Analytical Lab  
Page 7 of 28

*This report shall not be reproduced, except in full.*

Order # J11120115

Site: Trip Blank

Collection Date: 14-Dec-11 11:30 AM

Sample #: 2011026347

Matrix: OTHER

Analyte	Result	Units	Qualifiers	RDL	Method	Analysis Date/Time	Analyst
<b><u>TOTAL RECOVERABLE METALS BY ICP</u></b>							
Boron (B)	< 0.050	mg/L		0.05	EPA 200.7	16-Dec-11 10:13	DJSULL1
<b><u>TOTAL RECOVERABLE METALS BY ICP-MS</u></b>							
Arsenic (As)	< 1.000	ug/L		1	EPA 200.8	21-Dec-11 09:40	MHH7131
Chromium (Cr)	< 1.000	ug/L		1	EPA 200.8	21-Dec-11 09:40	MHH7131
Copper (Cu)	< 1.000	ug/L		1	EPA 200.8	21-Dec-11 09:40	MHH7131
Nickel (Ni)	< 1.000	ug/L		1	EPA 200.8	21-Dec-11 09:40	MHH7131
Selenium (Se)	< 1.000	ug/L		1	EPA 200.8	21-Dec-11 09:40	MHH7131
Silver (Ag)	< 1.000	ug/L		1	EPA 200.8	21-Dec-11 09:40	MHH7131
Zinc (Zn)	< 2.00	ug/L		2	EPA 200.8	21-Dec-11 09:40	MHH7131

## **SELENIUM SPECIATION**

Vendor Parameter Complete V\_AS&C

Site: BIOREACTOR 1 INF.

Collection Date: 14-Dec-11 11:40 AM

Sample #: 2011026356

Matrix: OTHER

Analyte	Result	Units	Qualifiers	RDL	Method	Analysis Date/Time	Analyst
<b><u>MERCURY 1631</u></b>							
Vendor Parameter	Complete				V_BRAND		

Site: HG BLANK BIOREACTOR 1 INF.

Collection Date: 14-Dec-11 11:40 AM

Sample #: 2011026357

Matrix: OTHER

Analyte	Result	Units	Qualifiers	RDL	Method	Analysis Date/Time	Analyst
<b><u>MERCURY 1631</u></b>							
Vendor Parameter	Complete				V_BRAND		

Site: BIOREACTOR 2 INF.

Collection Date: 14-Dec-11 11:50 AM

Sample #: 2011026358

Matrix: OTHER

Analyte	Result	Units	Qualifiers	RDL	Method	Analysis Date/Time	Analyst
<b><u>MERCURY 1631</u></b>							
Vendor Parameter	Complete				V_BRAND		

# Certificate of Laboratory Analysis

Analytical Lab  
Page 8 of 28

*This report shall not be reproduced, except in full.*

**Order # J11120115**

Site: Hg Blk BioReactor 2 Inf

Collection Date: 14-Dec-11 11:50 AM

**Sample #:** 2011026359

Matrix: OTHER

Analyte	Result	Units	Qualifiers	RDL	Method	Analysis Date/Time	Analyst
<b><u>MERCURY 1631</u></b>							
Vendor Parameter	Complete				V_BRAND		

Site: BIOREACTOR 2 EFF.

Collection Date: 14-Dec-11 11:45 AM

**Sample #:** 2011026360

Matrix: OTHER

Analyte	Result	Units	Qualifiers	RDL	Method	Analysis Date/Time	Analyst
<b><u>MERCURY 1631</u></b>							
Vendor Parameter	Complete				V_BRAND		

Site: Hg Blk BioReactor 2 Eff

Collection Date: 14-Dec-11 11:45 AM

**Sample #:** 2011026361

Matrix: OTHER

Analyte	Result	Units	Qualifiers	RDL	Method	Analysis Date/Time	Analyst
<b><u>MERCURY 1631</u></b>							
Vendor Parameter	Complete				V_BRAND		





**APPLIED SPECIATION  
AND CONSULTING, LLC**

18804 Northcreek Parkway Bothell, WA, 98011  
Tel: (425) 483-3300 Fax: (425) 483-9818  
[www.appliedspeciation.com](http://www.appliedspeciation.com)

December 22, 2011

Jay Perkins  
Duke Energy Analytical Laboratory  
Mail Code MGO3A2 (Building 7405)  
13339 Hagers Ferry Rd.  
Huntersville, NC 28078  
(704) 875-5245

Project: Belews – FGD WWTS Bi-Monthly Sampling) (LIMS # J11120115)

Dear Mr. Perkins,

Attached is the report associated with four (4) aqueous samples submitted for selenium speciation analysis on December 15, 2011. The samples were received in a sealed cooler at -0.7°C on December 16, 2011. Selenium speciation analysis was performed via ion chromatography inductively coupled plasma dynamic reaction cell mass spectrometry (IC-ICP-DRC-MS). Any issues associated with the analysis are addressed in the following report.

If you have any questions, please feel free to contact me at your convenience.

Sincerely,

A handwritten signature in black ink, appearing to read "Russell Gerads", with a large, stylized flourish at the end.

Russell Gerads  
Vice President  
Applied Speciation and Consulting, LLC

Applied Speciation and Consulting, LLC

Report prepared for:

Jay Perkins  
Duke Energy Analytical Laboratory  
Mail Code MGO3A2 (Building 7405)  
13339 Hagers Ferry Rd.  
Huntersville, NC 28078

Project: Belews – FGD WWTS Bi-Monthly Sampling) (LIMS # J11120115)

December 22, 2011

## 1. Sample Reception

Four (4) aqueous samples in 125mL HDPE bottles (provided by Applied Speciation and Consulting) were submitted for selenium speciation analysis on December 15, 2011. The samples were received on December 16, 2011 in a sealed container at -0.7°C.

The samples were received in a laminar flow clean hood, void of trace metals contamination and ultra-violet radiation, and was designated a discrete sample identifier. An aliquot of each sample was filtered (0.45µm) and each filtrate was stored in a secure, monitored cryofreezer (maintained at a temperature of -80°C) until selenium speciation analysis could be performed via ion chromatography inductively coupled plasma dynamic reaction cell mass spectrometry (IC-ICP-DRC-MS).

## 2. Sample Preparation

All sample preparation is performed in laminar flow clean hoods known to be free from trace metals contamination. All applied water for dilutions and sample preservatives are monitored for contamination to account for any biases associated with the sample results.

Selenium Speciation Analysis by IC-ICP-DRC-MS Prior to analysis, an aliquot of each sample was filtered with a syringe filter (0.45µm) and injected directly into a sealed autosampler vial. No further sample preparation was performed as any chemical alteration of a sample may shift the equilibrium of the system, resulting in changes in speciation ratios.

## 3. Sample Analysis

All sample analysis is preceded by a minimum of a five-point calibration curve spanning the entire concentration range of interest. Calibration curves are performed at the beginning of

each analytical day. All calibration curves, associated with each species of interest, are standardized by linear regression resulting in a response factor. All sample results are **instrument blank corrected** to account for any operational biases associated with the analytical platform.

Prior to sample analysis, all calibration curves are verified using second source standards which are identified as initial calibration verification standards (ICV).

Ongoing instrument performance is identified by the analysis of continuing calibration verification standards (CCV) and continuing calibration blanks (CCB) at a minimal interval of every ten analytical runs.

*Selenium Speciation Analysis by IC-ICP-DRC-MS* Each sample for selenium speciation analysis was analyzed by ion chromatography inductively coupled plasma dynamic reaction cell mass spectrometry (IC-ICP-DRC-MS) on December 16, 2011. An aliquot of each sample is injected onto an anion exchange column and mobilized by a basic ( $\text{pH} > 7$ ) gradient. The eluting selenium species are then introduced into a radio frequency (RF) plasma where energy-transfer processes cause desolvation, atomization, and ionization. The ions are extracted from the plasma through a differentially-pumped vacuum interface and travel through a pressurized chamber (DRC) containing a specific reactive gas which preferentially reacts with interfering ions of the same target mass to charge ratios ( $m/z$ ). A solid-state detector detects ions transmitted through the mass analyzer and the resulting current is processed by a data handling system.

Retention times for each eluting species are compared to known standards for species identification.

#### 4. Analytical Issues

The overall analyses went well and no significant analytical issues were encountered. All quality control parameters associated with this sample were within acceptance limits.

The estimated method detection limits (eMDLs) for selenite, selenate, and selenocyanate are generated from replicate analyses of the lowest standard in the calibration curve. Not all selenium species are present in preparation blanks; therefore, eMDL calculations based on preparation blanks are artificially biased low.

The eMDL for methylseleninic acid and selenomethionine is calculated from the average eMDL of selenite, selenate, and selenocyanate. The calibration does not contain methylseleninic acid or selenomethionine due to impurities in these standards which would bias the results for other selenium species.

If you have any questions or concerns regarding this report, please feel free to contact me.

Sincerely,

A handwritten signature in black ink, appearing to read 'Russell Gerads', with a stylized, flowing script.

Russell Gerads  
Vice President  
Applied Speciation and Consulting, LLC

Selenium Speciation Results for Duke Energy  
Project Name: Belews - FGD WWTS Bi-Monthly Sampling)  
Contact: Jay Perkins  
LIMS #J11120115

Date: December 22, 2011  
Report Generated by: Russell Gerads  
Applied Speciation and Consulting, LLC

**Sample Results**

Sample ID	Se(IV)	Se(VI)	SeCN	MeSe(IV)	SeMe	Unknown Se Species (n)
FGD Purge Eff	54.1	85.8	7.6	ND (<1.9)	ND (<1.9)	0 (0)
BioReactor 1 Inf	19.3	92.0	ND (<0.39)	ND (<0.47)	ND (<0.47)	2.00 (1)
BioReactor 2 Eff	ND (<0.56)	ND (<0.47)	ND (<0.39)	ND (<0.47)	ND (<0.47)	0 (0)
Metals Trip Blk	ND (<0.11)	ND (<0.093)	ND (<0.078)	ND (<0.094)	ND (<0.094)	0 (0)

All results reflect the applied dilution and are reported in µg/L

ND = Not detected at the applied dilution

SeCN = Selenocyanate

MeSe(IV) = Methylseleninic acid

SeMe = Selenomethionine

Unknown Se Species = Total concentration of all unknown Se species observed by IC-ICP-MS

n = number of unknown Se species observed

Selenium Speciation Results for Duke Energy  
Project Name: Belews - FGD WWTS Bi-Monthly Sampling)  
Contact: Jay Perkins  
LIMS #J11120115

Date: December 22, 2011  
Report Generated by: Russell Gerads  
Applied Speciation and Consulting, LLC

**Quality Control Summary - Preparation Blank Summary**

Analyte (µg/L)	PBW1	PBW2	PBW3	PBW4	Mean	StdDev	eMDL*	eMDL 10x	eMDL 50x	eMDL 200x
Se(IV)	0.000	0.000	0.000	0.000	0.000	0.000	0.011	0.11	0.56	2.2
Se(VI)	0.000	0.000	0.000	0.000	0.000	0.000	0.009	0.093	0.47	1.9
SeCN	0.000	0.000	0.000	0.000	0.000	0.000	0.008	0.078	0.39	1.6
MeSe(IV)	0.000	0.000	0.000	0.000	0.000	0.000	0.009	0.094	0.47	1.9
SeMe	0.000	0.000	0.000	0.000	0.000	0.000	0.009	0.094	0.47	1.9

eMDL = Estimated Method Detection Limit

\*Please see narrative regarding eMDL calculations

**Quality Control Summary - Certified Reference Materials**

Analyte (µg/L)	CRM	True Value	Result	Recovery
Se(IV)	LCS	3.83	4.09	106.8
Se(VI)	LCS	3.79	4.18	110.1
SeCN	LCS	3.57	3.83	107.3
MeSe(IV)	LCS	2.59	2.52	97.4
SeMe	LCS	3.73	3.77	101.1

Selenium Speciation Results for Duke Energy  
Project Name: Belews - FGD WWTS Bi-Monthly Sampling)  
Contact: Jay Perkins  
LIMS #J11120115

Date: December 22, 2011  
Report Generated by: Russell Gerads  
Applied Speciation and Consulting, LLC

**Quality Control Summary - Matrix Duplicates**

Analyte (µg/L)	Sample ID	Rep 1	Rep 2	Mean	RPD
Se(IV)	BioReactor 2 Eff	ND (<0.56)	ND (<0.56)	NC	NC
Se(VI)	BioReactor 2 Eff	ND (<0.47)	ND (<0.47)	NC	NC
SeCN	BioReactor 2 Eff	ND (<0.39)	ND (<0.39)	NC	NC
MeSe(IV)	BioReactor 2 Eff	ND (<0.47)	ND (<0.47)	NC	NC
SeMe	BioReactor 2 Eff	ND (<0.47)	ND (<0.47)	NC	NC

ND = Not detected at the applied dilution

NC = Value was not calculated due to one or more concentrations below the eMDL

**Quality Control Summary - Matrix Spike/ Matrix Spike Duplicate**

Analyte (µg/L)	Sample ID	Spike Conc	MS Result	Recovery	Spike Conc	MSD Result	Recovery	RPD
Se(IV)	BioReactor 2 Eff	278.0	266.0	95.7	278.0	261.6	94.1	1.6
Se(VI)	BioReactor 2 Eff	252.3	252.5	100.1	252.3	251.7	99.8	0.3
SeCN	BioReactor 2 Eff	228.8	211.8	92.6	228.8	213.0	93.1	0.6



Duke Energy Analytical Laboratory

Mail Code MGO3A2 (Building 7405)  
13338 Hagers Ferry Rd  
Huntersville, N.C. 28078  
(704) 875-5245  
Fax: (704) 875-4349

# CHAIN OF CUSTODY RECORD AND ANALYSIS REQUEST FORM

Customer must Complete

1) Project Name	Belwets - FGD		2) Phone No:
3) Client	WWTS Bi-Monthly Sampling)		4) Fax No:
5) Business Unit	Bill Kennedy, Melonie Martin, Wayne Chapman, Tom Johnson		6) Process:
8) Oper. Unit:	9) Res. Type:	10) Reso. Center:	

LAB USE ONLY
11) Lab ID
201026219
342
343
344
345
346
347

Customer to complete appropriate columns to right

Se Speciation Bottle ID	13 Sample Description or ID	Date	Time	Signature	17 Comp.	18 Grab	TDS	Hg - 245.1	Br (Dionex)	Metals*	Se, soluble (no dig.)	Se, speciation - vendor to AS&C (Important to place filled bottle back into both baggies)
	FGD Purge Eff	12/14	11:30	Thomas Thoms			1	1	1	1	1	1
	EQ Tank Eff.	12/14	11:30	Thomas Thoms			1	1	1	1	1	1
	BioReactor 1 Inf	12/14	11:30	Thomas Thoms								
	BioReactor 2 Inf	12/14	11:30	Thomas Thoms								
	BioReactor 2 Eff	12/14	11:30	Thomas Thoms								
	Filter Blk	12/14	11:30	Thomas Thoms								
	Metals Trip Blk	12/14	11:30	Thomas Thoms								

Filtering of the Se is performed in the field please provide a filter blank too.

Customer to sign & date below - fill out from left to right

1) Relinquished By	12/14	11:40	12/14/11	11:50
2) Relinquished By	12/14	11:50	12/14/11	12:25
3) Relinquished By	12/14	11:50	12/14/11	12:25
4) Relinquished By	12/14	11:50	12/14/11	12:25
5) Relinquished By	12/14	11:50	12/14/11	12:25
6) Relinquished By	12/14	11:50	12/14/11	12:25
7) Relinquished By	12/14	11:50	12/14/11	12:25
8) Seal/Locked By	12/14	11:50	12/14/11	12:25
9) Seal/Locked By	12/14	11:50	12/14/11	12:25
10) Seal/Locked By	12/14	11:50	12/14/11	12:25
11) Seal/Locked By	12/14	11:50	12/14/11	12:25

Customer, IMPORTANT!  
Please indicate desired turnaround.

22 Requested Turnaround

14 Days \_\_\_\_\_

7 Days \_\_\_\_\_

48 Hr \_\_\_\_\_

Other \_\_\_\_\_

\* Add Cost Will Apply

12-22-11

ORDER #	AS&C	PO#133241
DATE	12/14/11	11:50
LOGGED BY	Thomas Thoms	
DATE & TIME	12/14/11	11:50
MATRIX	OTHER	
COOLER TEMP (C)	4.34	
2-H <sub>2</sub> SO <sub>4</sub> 3-HNO <sub>3</sub>	4.34	
4-HCl 5-HNO <sub>3</sub>	4.34	
SAMPLE PROGRAM	Ground	
DRINKING WATER	Drinking Water	
ROTA Waste	ROTA Waste	

19 Page 1 of 2  
DISTRIBUTION  
ORIGINAL TO LAB,  
COPY TO CLIENT





December 28, 2011

Duke Energy  
ATTN: Jay Perkins  
13339 Hagers Ferry Road  
Huntersville NC 28078  
jperkins@duke-energy.com

RE: Project DUK-HV1101

Client Project: J11120115

Dear Jay Perkins,

This report contains results for the 6 samples received by Brooks Rand Labs (BRL) on December 16, 2011. The samples were logged-in for the contracted analyses according to the chain-of-custody form(s). The samples were received, prepared, analyzed, and stored according to BRL SOPs and EPA methodology.

The results were method blank corrected as described in the calculations section of the relevant BRL SOP(s) and may have been evaluated using reporting limits that have been adjusted to account for sample aliquot size. Please refer to the *Sample Results* page for sample-specific MDLs, MRLs, and other details. All data is reported without qualification (with the exception of concentration qualifiers), and all associated quality control sample results meet the acceptance criteria.

BRL, an accredited laboratory, certifies that the reported results of all analyses for which BRL is NELAP accredited meet all NELAP requirements. For more details, please see the *Report Information* page in your report. Please feel free to contact me if you have any questions regarding this report.

Sincerely,

Tiffany Stilwater  
Project Manager  
tiffany@brooksrands.com

## Report Information

### Laboratory Accreditation

BRL is accredited by the *National Environmental Laboratory Accreditation Program* (NELAP) through the State of Florida Department of Health, Bureau of Laboratories (E87982) and is certified to perform many environmental analyses. BRL is also certified by many other states to perform environmental analyses. For a current list of our accreditations/certifications, please visit our website at <http://www.brooksrand.com/default.asp?contentID=586>. Results reported relate only to the samples listed in the report.

### Field Quality Control Samples

Please be notified that certain EPA methods require the collection of field quality control samples of an appropriate type and frequency; failure to do so is considered a deviation from some methods and for compliance purposes should only be done with the approval of regulatory authorities. Please see the specific EPA methods for details regarding required field quality control samples.

### Common Abbreviations

<b>BLK</b>	method blank	<b>MS</b>	matrix spike
<b>BRL</b>	Brooks Rand Labs	<b>MSD</b>	matrix spike duplicate
<b>BS</b>	laboratory fortified blank	<b>ND</b>	non-detect
<b>CAL</b>	calibration standard	<b>NR</b>	non-reportable
<b>CCV</b>	continuing calibration verification	<b>PS</b>	post preparation spike
<b>COC</b>	chain of custody record	<b>REC</b>	percent recovery
<b>CRM</b>	certified reference material	<b>RPD</b>	relative percent difference
<b>D</b>	dissolved fraction	<b>RSD</b>	relative standard deviation
<b>DUP</b>	duplicate	<b>SCV</b>	secondary calibration verification
<b>ICV</b>	initial calibration verification	<b>SOP</b>	standard operating procedure
<b>MDL</b>	method detection limit	<b>SRM</b>	standard reference material
<b>MRL</b>	method reporting limit	<b>T</b>	total recoverable fraction

### Definition of Data Qualifiers

(Effective 9/23/09)

<b>B</b>	Detected by the instrument, the result is > the MDL but ≤ the MRL. Result is reported and considered an estimate.
<b>E</b>	An estimated value due to the presence of interferences. A full explanation is presented in the narrative.
<b>H</b>	Holding time and/or preservation requirements not met. Result is estimated.
<b>J</b>	Estimated value. A full explanation is presented in the narrative.
<b>J-M</b>	Duplicate precision (RPD) for associated QC sample was not within acceptance criteria. Result is estimated.
<b>J-N</b>	Spike recovery for associated QC sample was not within acceptance criteria. Result is estimated.
<b>M</b>	Duplicate precision (RPD) was not within acceptance criteria. Result is estimated.
<b>N</b>	Spike recovery was not within acceptance criteria. Result is estimated.
<b>R</b>	Rejected, unusable value. A full explanation is presented in the narrative.
<b>U</b>	Result is ≤ the MDL or client requested reporting limit (CRRL). Result reported as the MDL or CRRL.
<b>X</b>	Result is not BLK-corrected and is within 10x the absolute value of the highest detectable BLK in the batch. Result is estimated.

These qualifiers are based on those previously utilized by Brooks Rand, Ltd., those found in the EPA SOW ILM03.0, Exhibit B, Section III, pg. B-18, and the USEPA Laboratory Data Validation Functional Guidelines for Evaluating Inorganic Analyses; USEPA; July 2002. These supersede all previous qualifiers ever employed by BRL.

**Project ID:** DUK-HV1101  
**PM:** Tiffany Stilwater



Analytical Lab  
Page 19 of 28  
Report 1151039  
**Client PM:** Jay Perkins  
**Client PO:** 141391

## Sample Information

Sample	Lab ID	Report Matrix	Type	Sampled	Received
BioReactor 1 Inf	1151039-01	Influent	Sample	12/14/2011	12/16/2011
Hg Blk BioReactor 1 Inf	1151039-02	DIW	Field Blank	12/14/2011	12/16/2011
BioReactor 2 Inf	1151039-03	Influent	QC Sample	12/14/2011	12/16/2011
Hg Blk BioReactor 2 Inf	1151039-04	DIW	Field Blank	12/14/2011	12/16/2011
BioReactor 2 Eff	1151039-05	Effluent	Sample	12/14/2011	12/16/2011
Hg Blk BioReactor 2 Eff	1151039-06	DIW	Field Blank	12/14/2011	12/16/2011

## Batch Summary

Analyte	Lab Matrix	Method	Prepared	Analyzed	Batch	Sequence
Hg	Water	EPA 1631	12/19/2011	12/20/2011	B112132	1100906



## Sample Results

Sample	Analyte	Report Matrix	Fraction	Result	Qualifier	MDL	MRL	Unit	Batch	Sequence
<b>BioReactor 1 Inf</b>										
1151039-01	Hg	Influent	T	161		1.52	4.04	ng/L	B112132	1100906
<b>BioReactor 2 Eff</b>										
1151039-05	Hg	Effluent	T	39.4		0.61	1.61	ng/L	B112132	1100906
<b>BioReactor 2 Inf</b>										
1151039-03	Hg	Influent	T	102		1.52	4.04	ng/L	B112132	1100906
<b>Hg Blk BioReactor 1 Inf</b>										
1151039-02	Hg	DIW	T	0.15	U	0.15	0.40	ng/L	B112132	1100906
<b>Hg Blk BioReactor 2 Eff</b>										
1151039-06	Hg	DIW	T	0.15	U	0.15	0.40	ng/L	B112132	1100906
<b>Hg Blk BioReactor 2 Inf</b>										
1151039-04	Hg	DIW	T	0.15	U	0.15	0.40	ng/L	B112132	1100906

**Project ID:** DUK-HV1101  
**PM:** Tiffany Stilwater



Analytical Lab  
 Page 2 of 28  
 Report 1151039  
**Client PM:** Jay Perkins  
**Client PO:** 141391

## Accuracy & Precision Summary

**Batch:** B112132  
**Lab Matrix:** Water  
**Method:** EPA 1631

Sample	Analyte	Native	Spike	Result	Units	REC & Limits	RPD & Limits
<b>B112132-SRM1</b>	<b>Certified Reference Material (1149037, NIST 1641d 1000x dilution)</b>						
	Hg		15.68	15.29	ng/L	98% 85-115	
<b>B112132-MS2</b>	<b>Matrix Spike (1151039-03)</b>						
	Hg	102.3	505.1	664.3	ng/L	111% 71-125	
<b>B112132-MSD2</b>	<b>Matrix Spike Duplicate (1151039-03)</b>						
	Hg	102.3	505.1	648.2	ng/L	108% 71-125	2% 24
<b>B112132-MS3</b>	<b>Matrix Spike (1151039-05)</b>						
	Hg	39.35	177.4	228.6	ng/L	107% 71-125	
<b>B112132-MSD3</b>	<b>Matrix Spike Duplicate (1151039-05)</b>						
	Hg	39.35	180.9	226.1	ng/L	103% 71-125	1% 24

**Project ID:** DUK-HV1101  
**PM:** Tiffany Stilwater



Analytical Lab  
Page 22 of 28  
Report 1151039  
**Client PM:** Jay Perkins  
**Client PO:** 141391

## Method Blanks & Reporting Limits

**Batch:** B112132  
**Matrix:** Water  
**Method:** EPA 1631  
**Analyte:** Hg

Sample	Result	Units
B112132-BLK1	0.06	ng/L
B112132-BLK2	0.08	ng/L
B112132-BLK3	0.06	ng/L
B112132-BLK4	0.04	ng/L

**Average:** 0.06  
**Limit:** 0.50

**Standard Deviation:** 0.02  
**Limit:** 0.10

**MDL:** 0.15  
**MRL:** 0.40

**Project ID:** DUK-HV1101  
**PM:** Tiffany Stilwater



Analytical Lab  
Page 23 of 28  
Report 1151039  
**Client PM:** Jay Perkins  
**Client PO:** 141391

## Instrument Calibration

**Sequence:** 1100906  
**Instrument:** THG-10  
**Date:** 12/20/2011  
**Analyte:** Hg

**Total Mercury and Mercury Speciation by CVAFS**  
**Method:** EPA 1631

Lab ID	True Value	Result	Units	REC & Limits	
1100906-IBL1		3.72	pg of Hg		
1100906-IBL2		3.49	pg of Hg		
1100906-IBL3		3.46	pg of Hg		
1100906-IBL4		4.98	pg of Hg		
1100906-CAL1	25.00	24.38	pg of Hg	98%	
1100906-CAL2	100.0	98.22	pg of Hg	98%	
1100906-CAL3	500.0	482.6	pg of Hg	97%	
1100906-CAL4	2500	2687	pg of Hg	107%	
1100906-CAL5	10000	10100	pg of Hg	101%	
1100906-ICV1	1568	1529	pg of Hg	98%	85-115
1100906-CCB1		6.98	pg of Hg		
1100906-CCV1	500.0	532.2	pg of Hg	106%	77-123
1100906-CCV2	500.0	515.4	pg of Hg	103%	77-123
1100906-CCV3	500.0	514.5	pg of Hg	103%	77-123

**Project ID:** DUK-HV1101  
**PM:** Tiffany Stilwater



Analytical Lab  
Page 24 of 28  
Report 1151039  
**Client PM:** Jay Perkins  
**Client PO:** 141391

## Sample Containers

<b>Lab ID:</b> 1151039-01			<b>Report Matrix:</b> Influent			<b>Collected:</b> 12/14/2011	
<b>Sample:</b> BioReactor 1 Inf			<b>Sample Type:</b> Sample			<b>Received:</b> 12/16/2011	
<b>Des</b>	<b>Container</b>	<b>Size</b>	<b>Lot</b>	<b>Preservation</b>	<b>P-Lot</b>	<b>pH</b>	<b>Ship. Cont.</b>
A	Bottle FLPE Hg-T	250 mL	71470160 10	none	n/a		Cooler
<b>Lab ID:</b> 1151039-02			<b>Report Matrix:</b> DIW			<b>Collected:</b> 12/14/2011	
<b>Sample:</b> Hg Blk BioReactor 1 Inf			<b>Sample Type:</b> Field Blank			<b>Received:</b> 12/16/2011	
<b>Des</b>	<b>Container</b>	<b>Size</b>	<b>Lot</b>	<b>Preservation</b>	<b>P-Lot</b>	<b>pH</b>	<b>Ship. Cont.</b>
A	Bottle FLPE Hg-T	250 mL	71470160 10	none	n/a		Cooler
<b>Lab ID:</b> 1151039-03			<b>Report Matrix:</b> Influent			<b>Collected:</b> 12/14/2011	
<b>Sample:</b> BioReactor 2 Inf			<b>Sample Type:</b> QC Sample			<b>Received:</b> 12/16/2011	
<b>Des</b>	<b>Container</b>	<b>Size</b>	<b>Lot</b>	<b>Preservation</b>	<b>P-Lot</b>	<b>pH</b>	<b>Ship. Cont.</b>
A	Bottle FLPE Hg-T	250 mL	71470160 10	none	n/a		Cooler
<b>Lab ID:</b> 1151039-04			<b>Report Matrix:</b> DIW			<b>Collected:</b> 12/14/2011	
<b>Sample:</b> Hg Blk BioReactor 2 Inf			<b>Sample Type:</b> Field Blank			<b>Received:</b> 12/16/2011	
<b>Des</b>	<b>Container</b>	<b>Size</b>	<b>Lot</b>	<b>Preservation</b>	<b>P-Lot</b>	<b>pH</b>	<b>Ship. Cont.</b>
A	Bottle FLPE Hg-T	250 mL	71470160 10	none	n/a		Cooler
<b>Lab ID:</b> 1151039-05			<b>Report Matrix:</b> Effluent			<b>Collected:</b> 12/14/2011	
<b>Sample:</b> BioReactor 2 Eff			<b>Sample Type:</b> Sample			<b>Received:</b> 12/16/2011	
<b>Des</b>	<b>Container</b>	<b>Size</b>	<b>Lot</b>	<b>Preservation</b>	<b>P-Lot</b>	<b>pH</b>	<b>Ship. Cont.</b>
A	Bottle FLPE Hg-T	250 mL	71470160 10	none	n/a		Cooler
<b>Lab ID:</b> 1151039-06			<b>Report Matrix:</b> DIW			<b>Collected:</b> 12/14/2011	
<b>Sample:</b> Hg Blk BioReactor 2 Eff			<b>Sample Type:</b> Field Blank			<b>Received:</b> 12/16/2011	
<b>Des</b>	<b>Container</b>	<b>Size</b>	<b>Lot</b>	<b>Preservation</b>	<b>P-Lot</b>	<b>pH</b>	<b>Ship. Cont.</b>
A	Bottle FLPE Hg-T	250 mL	71470160 10	none	n/a		Cooler



**Project ID:** DUK-HV1101  
**PM:** Tiffany Stilwater



Analytical Lab  
Page 25 of 28  
Report 1151039  
**Client PM:** Jay Perkins  
**Client PO:** 141391

## Shipping Containers

### Cooler

**Received:** December 16, 2011 9:00  
**Tracking No:** 4726 7966 6868 via FedEx  
**Coolant Type:** Ice  
**Temperature:** 1.7 °C

**Description:** Cooler  
**Damaged in transit?** No  
**Returned to client?** No

**Custody seals present?** No  
**Custody seals intact?** No  
**COC present?** Yes

# CHAIN OF CUSTODY RECORD AND ANALYSIS REQUEST FORM

1151039  
Analytical Lab  
Page 2 of 2  
1151039



**Duke Energy Analytical Laboratory**  
Mail Code MGO3A2 (Building 7405)  
13339 Hagers Ferry Rd  
Huntersville, N. C. 28078  
(704) 875-5245  
Fax: (704) 875-4349

Analytical Laboratory Use Only			
ORDER # J 11/20/15	Sample Class OTHER	Samples Originating From	NC SC
Logged By BA	Date & Time 12/14/11 1536	SAMPLE PROGRAM Water _____ Ground NPDES Drinking Water UST RCRA Waste _____	
Vt	Brooks Rand PO#141391	Cooler Temp (C) Preserv.: 1=HCL 2=H <sub>2</sub> SO <sub>4</sub> 3=HNO <sub>3</sub> 4=Ice 5=None	

19 Page 2 of 2  
DISTRIBUTION  
ORIGINAL to LAB,  
COPY to CLIENT

Customer must Complete

1) Project Name <b>Belews - FGD</b>	2) Phone No:
2) Client: <b>WWTS (2011, Bi-Weekly Sampling)</b>	4) Fax No:
5) Business Unit:	6) Process:
8) Oper. Unit:	10) Reso. Center:
Mail Code:	

Customer to complete all appropriate non-shaded areas.						16 Analyses Required		5	
Sampling conducted: 2nd Wednesday each month						17 Comp.	18 Grab	Hg 1637 (sample 2nd week only)	
13 Sample Description or ID	Date	Time	Signature						
BioReactor 1 Inf	12-14-11	1140	<i>David M...</i>						
Hg Blk BioReactor 1 Inf								1	
BioReactor 2 Inf		1150						1	
Hg Blk BioReactor 2 Inf								1	
BioReactor 2 Eff		1145						1	
Hg Blk BioReactor 2 Eff								1	

Use the Bioreactor 2 Inf or EFF sample as the MS/MSD

Customer to complete appropriate columns to right

LAB USE ONLY
11 Lab ID
2011026356
57
58
59
60
61

1) Relinquished By: <i>David M...</i>	Date/Time: 12-14-11 1525	2) Accepted By: <i>Cal M...</i>	Date/Time: 12-14-11 1525
3) Relinquished By:	Date/Time:	4) Accepted By: <i>[Signature]</i>	Date/Time: 12/16/11 0900
5) Relinquished By:	Date/Time:	6) Accepted By:	Date/Time:
7) Relinquished By:	Date/Time:	8) Accepted By:	Date/Time:
9) Seal/Locked By:	Date/Time:	10) Seal/Lock Opened By:	Date/Time:
11) Seal/Locked By:	Date/Time:	12) Seal/Lock Opened By:	Date/Time:
Comments			

\* Metals=As, Ag, B, Cu, Cr, Ni, Se, Zn \*thomas.d.johnson@siemens.com

Customer, IMPORTANT!  
Please indicate desired turnaround.

22 Requested Turnaround
14 Days _____
*7 Days _____
*48 Hr _____
*Other _____
* Add. Cost Will Apply
12-22-11





## Duke Energy Analytical Laboratory

**Mail Code MGO3A2 (Building 7405)**  
**13339 Hagers Ferry Rd**  
**Huntersville, N. C. 28078**  
**(704) 875-5245**  
**Fax: (704) 875-4349**

## Analytical Laboratory Use Only

ORDER# J11120115	MATRIX: OTHER	Samples Originating From	NC SC
Logged By RA	Date & Time 12/14/11 1536	SAMPLE PROGRAM Water	Ground NPDES Drinking Water UST RCRA Waste
AS&C	21.0 Cooler Temp (C)		

19 Page 1 of 2 Analytical Lab  
Page 27 of 28  
**DISTRIBUTION**  
ORIGINAL to LAB,  
COPY to CLIENT

1)Project Name		2)Phone No:
<b>Belews - FGD</b> <b>WWTS Bi-Monthly Sampling)</b>		
2) Client:		4)Fax No:
<b>Bill Kennedy, Melonie Martin,</b> <b>Wayne Chapman, Tom Johnson **</b>		
5)Business Unit:	6)Process:	Mail Code:
8)Oper. Unit:	9)Res. Type:	10)Reso. Center:

PO#133241			<sup>15</sup> Preserv.: 1=HCL 2=H <sub>2</sub> SO <sub>4</sub> 3=HNO <sub>3</sub> 4=Ice 5=None		4	3,4			4	3,4						4
MR #			<sup>16</sup> Analyses Required													
Customer to complete all appropriate non-shaded areas.																
Sampling conducted: 2nd and 4th Wednesday																
Date	Time	Signature	<sup>17</sup> Comp.	<sup>18</sup> Grab	TDS	Hg - 245.1	Br (Dionex)	Metals*	Se, soluble (no dig.)							Se, speciation - vendor to AS&C (important to place filled bottle back into both baggies)
12/14	11:30	Travis Thoms			1	1	1	1	1							1
12/14	11:30	Travis Thoms				1		1	1							
12/14	11:30	Travis Thoms						1	1							1
12/14	11:30	Travis Thoms						1								
12/14	11:30	Travis Thoms				1	1	1								1
12/14	11:30	Travis Thoms							1							1
Filtering of the Se is performed in the field please provide a filter blank too.																

LAB USE ONLY

<sup>11</sup>Lab ID

2011026219

342

343

344

345

346

347

Complete appropriate columns to right

Customer to sign & date below - fill out from left to right.

Customer to sign & date below - Fill out from left to right			
1) Relinquished By <i>Tom Johnson</i>	Date/Time 12/14 11:40	2) Accepted By <i>David Morris</i>	Date/Time 12-14-11 11:40
3) Relinquished By <i>David Morris</i>	Date/Time 12-14-11 1525	4) Accepted By <i>GC Sharma</i>	Date/Time 12-14-11 1525
5) Relinquished By <i>R. Harris</i>	Date/Time 12/15/11 1300	6) Accepted By:	Date/Time
7) Relinquished By	Date/Time	8) Accepted By:	Date/Time
9) Seal/Locked By <i>R. Harris</i>	Date/Time 12/15/11 1300	10) Seal/Lock Opened By	Date/Time
11) Seal/Locked By	Date/Time	12) Seal/Lock Opened By	Date/Time
Comments			
* B by ICP      As, Cr, Cu, Ni, Se, Ag, Zn by IMS      Digestions = TRM <a href="mailto:thomas.d.johnson@siemens.com">thomas.d.johnson@siemens.com</a>			

**Customer, IMPORTANT!**  
Please indicate desired turnaround.

**22 Requested Turnaround**

14 Days \_\_\_\_\_

\*7 Days \_\_\_\_\_


\* 48 Hr \_\_\_\_\_

\*Other \_\_\_\_\_

\* Add. Cost Will Apply



Analytical Lab



Analytical Lab  
19 **Page 28 of 28**  
**DISTRIBUTION**  
ORIGINAL to LAB  
COPY to CLIENT

[illegible]

Customer to complete appropriate columns to right

**Customer, IMPORTANT!**  
Please indicate desired turnaround.

**22 Requested Turnaround**

14 Days \_\_\_\_\_

\*7 Days \_\_\_\_\_

\* 48 Hr \_\_\_\_\_

\*Other \_\_\_\_\_

\* Add. Cost Will Apply

\* Metals=As, Ag, B, Cu, Cr, Ni, Se, Zn      \*thomas.d.johnson@siemens.com