

checked by TT 2/24/17

CETIFICATION

SDG No:

JC36372

Laboratory:

Accutest, New Jersey

Site:

BMSMC, Humacao, PR

Matrix:

Soil

SUMMARY:

Soil samples (Table 1) were collected on the BMSMC facility. The BMSMC facility is located in Humacao, PR. Samples were taken January 26, 2017 and were analyzed in Accutest Laboratory of Dayton, New Jersey that reported the data under SDG No.: JC36372. Results were validated using the latest validation guidelines (July, 2015) of the EPA Hazardous Waste Support Section. The analyses performed are shown in Table 1. Individual data review worksheets are enclosed for each target analyte group. The data sample organic data samples summary form shows for analytes results that were qualified.

In summary the results are valid and can be used for decision taking purposes.

Table 1. Samples analyzed and analysis performed

SAMPLE ID	SAMPLE DESCRIPTION	MATRIX	ANALYSIS PERFORMED
JC36372-1	BKGSS-1	Soil	SVOCs; PAHs (SIM); LMWA; Pesticides; Metals
JC36372-2	BKGSS-2	Soil	SVOCs; PAHs (SIM); LMWA; Pesticides; Metals
JC36372-3	BKGSS-2 DUP	Soil	LMWA; Pesticides
JC36372-4	BKGSS-3	Soil	SVOCs; PAHs (SIM); LMWA; Pesticides; Metals
JC36372-4D	BKGSS-3 MSD	Soil	SVOCs; PAHs (SIM); LMWA; Pesticides; Metals
JC36372-4S	BKGSS-3 MS	Soil	SVOCs; PAHs (SIM); LMWA; Pesticides; Metals
JC36372-5	FB-012617	AQ – Field Blank Soil	SVOCs; PAHs (SIM); LMWA; Pesticides; Metals
JC36372-6	EB-012617	AQ – Equipment Blank	SVOCs; PAHs (SIM); LMWA; Pesticides; Metals

ifael Infat Méndez

1617180

Reviewer Name:

Rafael Infante

Chemist License 1888

Signature:

Date:

February 18, 2017

Report of Analysis

Page 1 of 3

Client Sample ID: BKG SS-1 Lab Sample ID:

JC36372-1

SO - Soil

SW846 8270D SW846 3546

Date Received:

Date Sampled: 01/26/17 01/27/17

Method: Project:

Matrix:

BMSMC, PR

Percent Solids: 86.2

File ID 6P34796.D Run #1

DF 1

By Analyzed 02/08/17 AC

Prep Date 01/30/17

Prep Batch OP173

Analytical Batch E6P1601

Run #2

Initial Weight

30.0 g

Final Volume 1.0 ml

Run #1 Run #2

ABN Special List

CAS No.	Compound	Result	RL	MDL	Units
95-57-8	2-Chlorophenol	ND	77	19	ug/kg
59-50-7	4-Chloro-3-methyl phenol	ND	190	24	ug/kg
120-83-2	2,4-Dichlorophenol	ND	190	33	ug/kg
105-67-9	2,4-Dimethylphenol	ND	190	69	ug/kg
51-28-5	2,4-Dinitrophenol	ND	190	150	ug/kg
534-52-1	4,6-Dinitro-o-cresol	ND	190	41	ug/kg
95-48-7	2-Methylphenol	ND	77	25	ug/kg
	3&4-Methylphenol	ND	77	32	ug/kg
88-75-5	2-Nitrophenol	ND	190	26	ug/kg
100-02-7	4-Nitrophenol	ND	390	100	ug/kg
87-86-5	Pentachlorophenol	ND	150	36	ug/kg
108-95-2	Phenol	ND	77	20	ug/kg
58-90-2	2,3,4,6-Tetrachlorophenol	ND	190	26	ug/kg
95-95-4	2,4,5-Trichlorophenol	ND	190	29	ug/kg
88-06-2	2,4,6-Trichlorophenol	ND	190	23	ug/kg
83-32-9	Acenaphthene	48.4	39	13	ug/kg
208-96-8	Acenaphthylene	ND	39	20	ug/kg
98-86-2	Acetophenone	ND	190	8.3	ug/kg
120-12-7	Anthracene	144	39	24	ug/kg
1912-24-9	Atrazine	ND	77	17	ug/kg
56-55-3	Benzo(a)anthracene	628	39	11	ug/kg
50-32-8	Benzo(a)pyrene	632	39	18	ug/kg
205-99-2	Benzo(b)fluoranthene	857	39	17	ug/kg
191-24-2	Benzo(g,h,i)perylene	557	39	19	ug/kg
207-08-9	Benzo(k)fluoranthene	278	39	18	ug/kg
101-55-3	4-Bromophenyl phenyl ether	ND	77	15	ug/kg
85-68-7	Butyl benzyl phthalate	ND	77	9.4	ug/kg
92-52-4	I, 1'-Biphenyl	ND	77	5.3	ug/kg
100-52-7	Benzaldehyde	31.6	190	9.6	ug/kg
91-58-7	2-Chloronaphthalene	ND	77	9.2	ug/kg
106-47-8	4-Chloroaniline	ND	190	14	ug/kg
86-74-8	Carbazole	85.9	77	5.6	ug/kg

ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

Page 2 of 3

Client Sample ID: BKG SS-1 Lab Sample ID:

JC36372-1

SO - Soil

Date Sampled: 01/26/17

Matrix: Method:

SW846 8270D SW846 3546

Percent Solids: 86.2

Q

J

Date Received: 01/27/17

Project:

BMSMC, PR

ABN Special List

CAS No.	Compound	Result	RL	MDL	Units
105-60-2	Caprolactam	ND	77	15	ug/kg
218-01-9	Chrysene	620	39	12	ug/kg
111-91-1	bis(2-Chloroethoxy)methane	ND	77	8.3	ug/kg
111-44-4	bis(2-Chloroethyl)ether	ND	77	17	ug/kg
108-60-1	bis(2-Chloroisopropyl)ether	ND	77	14	ug/kg
7005-72-3	4-Chlorophenyl phenyl ether	ND	77	13	ug/kg
121-14-2	2,4-Dinitrotoluene	ND	39	12	ug/kg
606-20-2	2,6-Dinitrotoluene	ND	39	19	ug/kg
91-94-1	3,3'-Dichlorobenzidine	ND	77	32	ug/kg
123-91-1	1,4-Dioxane	ND	39	26	ug/kg
132-64-9	Dibenzofuran	21.0	77	16	ug/kg
84-74-2	Di-n-butyl phthalate	ND	77	6.3	ug/kg
117-84-0	Di-n-octyl phthalate	ND	77	9.6	ug/kg
84-66-2	Diethyl phthalate	ND	77	8.2	ug/kg
131-11-3	Dimethyl phthalate	ND	77	6.9	ug/kg
117-81-7	bis(2-Ethylhexyl)phthalate	ND	77	9.0	ug/kg
206-44-0	Fluoranthene	1250	39	17	ug/kg
86-73-7	Fluorene	34.6	39	18	ug/kg
118-74-1	Hexachlorobenzene	ND	77	9.8	ug/kg
87-68-3	Hexachlorobutadiene	ND	39	16	ug/kg
77-47-4	Hexachlorocyclopentadiene	ND	390	15	ug/kg
67-72-1	Hexachloroethane	ND	190	19	ug/kg
193-39-5	Indeno(1,2,3-cd)pyrene	556	39	18	ug/kg
78-59-1	Isophorone	ND	77	8.3	ug/kg
90-12-0	I-Methylnaphthalene	ND	77	7.6	ug/kg
91-57-6	2-Methylnaphthalene	ND	77	8.7	ug/kg
88-74-4	2-Nitroaniline	ND	190	9.1	ug/kg
99-09-2	3-Nitroaniline	ND	190	9.7	ug/kg
100-01-6	4-Nitroaniline	ND	190	10	ug/kg
91-20-3	Naphthalene	129	39	11	ug/kg
98-95-3	Nitrobenzene	ND	77	15	ug/kg
621-64-7	N-Nitroso-di-n-propylamine	ND	77	11	ug/kg
86-30-6	N-Nitrosodiphenylamine	ND	190	14	ug/kg
85-01-8	Phenanthrene	661	39	13	ug/kg
129-00-0	Pyrene	1200	39	12	ug/kg
95-94-3	1,2,4,5-Tetrachlorobenzene	ND	190	9.8	ug/kg
CAS No	Surrogata Dagovarias	Dun# 1	D.,# 2	Limi	i de en

CAS No. **Surrogate Recoveries** Run# 1 Run# 2 Limits

2-Fluorophenol

367-12-4

MDL = Method Detection Limit

71%

RL = Reporting Limit

ND = Not detected

E = Indicates value exceeds calibration range

J = Indicates an estimated value

23-115%

B = Indicates analyte found in associated method blank

ifael Infante Méndez IC # 1888







Page 3 of 3

Client Sample ID: BKG SS-1 Lab Sample ID:

Matrix: Method:

Project:

JC36372-1 SO - Soil

BMSMC, PR

SW846 8270D SW846 3546


Date Sampled: 01/26/17 Date Received: 01/27/17

Percent Solids: 86.2

ABN Special List

CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limits
4165-62-2	Phenol-d5	69%		27-114%
118-79-6	2,4,6-Tribromophenol	77%		19-152%
4165-60-0	Nitrobenzene-d5	69%		26-134%
321-60-8	2-Fluorobiphenyl	71%		39-124%
1718-51-0	Terphenyl-d14	69%		36-134%





MDL = Method Detection Limit

J = Indicates an estimated value

RL = Reporting Limit

B = Indicates analyte found in associated method blank

E = Indicates value exceeds calibration range



Report of Analysis

Page 1 of 1

Client Sample ID: BKG SS-1 Lab Sample ID:

JC36372-1

Matrix:

SO - Soil

Dibenzo(a,h)anthracene

Method: Project:

SW846 8270D BY SIM SW846 3546

BMSMC, PR

Date Sampled: 01/26/17

Q

Date Received: 01/27/17

Percent Solids: 86.2

Prep Batch Analytical Batch File ID DF Analyzed By **Prep Date** Run #1 01/30/17 OP173A E4M3201 4M69657.D 1 01/30/17 SG

Run #2

53-70-3

4165-60-0

321-60-8

1718-51-0

Initial Weight 30.0 g

Final Volume 1.0 ml

Run #1 Run #2

Compound CAS No.

Result

RL

3.9

MDL

0.90

Units

ug/kg

CAS No. Surrogate Recoveries

Nitrobenzene-d5

2-Fluorobiphenyl

Terphenyl-d14

Run# 1

93.2

Run#2 Limits 15-138%

52% 53% 72%

12-148% 10-157%



ND = Not detected

MDL = Method Detection Limit

J = Indicates an estimated value

RL = Reporting Limit

B = Indicates analyte found in associated method blank

E = Indicates value exceeds calibration range

Report of Analysis

Page 1 of 1

Client Sample ID: BKG SS-1

Lab Sample ID: JC36372-1 Matrix: SO - Soil

Method:

SW846-8015C (DAI)

Project:

BMSMC, PR

Date Sampled: 01/26/17

Q

Date Received: 01/27/17

Percent Solids: 86.2

File ID DF Analyzed Ву **Prep Date Prep Batch Analytical Batch** XPL GGH5640 Run #1 GH108422.D 1 01/30/17 n/a n/a

Run #2

Initial Weight

Run #1 5.0 g

Run #2

Low Molecular Alcohol List

CAS No.	Compound	Result	RL	MDL	Units
	•				
64-17-5	Ethanol	ND	120	80	ug/kg
78-83-1	Isobutyl Alcohol	ND	120	68	ug/kg
67-63-0	Isopropyl Alcohol	ND	120	66	ug/kg
71-23-8	n-Propyl Alcohol	ND	120	47	ug/kg
71-36-3	n-Butyl Alcohol	ND	120	63	ug/kg
78-92-2	sec-Butyl Alcohol	ND	120	62	ug/kg
67-56-1	Methanol	ND	230	56	ug/kg

CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limits
111-27-3	Hexanol	95%		52-141%
111-27-3	Hexanol	93%		52-141%



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank



By

KD

Page 1 of 1

Client Sample ID: BKG SS-1

Lab Sample ID:

JC36372-1

Matrix: Method: SO - Soil

SW846 8081B SW846 3546

Analyzed

02/01/17

Date Sampled: 01/26/17 Date Received: 01/27/17

OP215

Q

Percent Solids: 86.2

Project:

BMSMC, PR

Prep Date

01/31/17

Prep Batch Analytical Batch

G4G2017

Run #1 Run #2

Run #1

Run #2

Initial Weight

Final Volume

15.3 g

File ID

4G77547.D

10.0 ml

DF

1

Pesticide TCL List

CAS No.	Compound	Result	RL	MDL	Units
309-00-2	Aldrin	ND	0.76	0.36	ug/kg
319-84-6	alpha-BHC	ND	0.76	0.41	ug/kg
319-85-7	beta-BHC	ND	0.76	0.48	ug/kg
319-86-8	delta-BHC	ND	0.76	0.34	ug/kg
58-89-9	gamma-BHC (Lindane)	ND	0.76	0.34	ug/kg
5103-71-9	alpha-Chlordane	ND	0.76	0.36	ug/kg
5103-74-2	gamma-Chlordane	ND	0.76	0.33	ug/kg
60-57-1	Dieldrin	ND	0.76	0.38	ug/kg
72-54-8	4,4'-DDD	ND	0.76	0.49	ug/kg
72-55-9	4,4'-DDE	ND	0.76	0.40	ug/kg
50-29-3	4,4'-DDT	ND	0.76	0.45	ug/kg
72-20-8	Endrin	ND	0.76	0.36	ug/kg
1031-07-8	Endosulfan sulfate	ND	0.76	0.30	ug/kg
7421-93-4	Endrin aldehyde	ND	0.76	0.45	ug/kg
959-98-8	Endosulfan-l	ND	0.76	0.40	ug/kg
33213-65-9	Endosulfan-II	ND	0.76	0.40	ug/kg
76-44-8	Heptachlor	ND	0.76	0.37	ug/kg
1024-57-3	Heptachlor epoxide	ND	0.76	0.41	ug/kg
72-43-5	Methoxychlor	ND	1.5	0.38	ug/kg
53494-70-5	Endrin ketone	ND	0.76	0.58	ug/kg
8001-35-2	Toxaphene	ND	19	7.9	ug/kg
CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Lim	its

877-09-8	Tetrachloro-m-xylene	84%	25-135%
877-09-8	Tetrachloro-m-xylene	78%	25-135%
2051-24-3	Decachlorobiphenyl	130%	10-156%
2051-24-3	Decachlorobiphenyl	86%	10-156%



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank



Page 1 of 1

Client Sample ID: BKG SS-1 Lab Sample ID: Matrix:

JC36372-1 SO - Soil

Date Sampled: 01/26/17 Date Received: 01/27/17

Percent Solids: 86.2

Project:

BMSMC, PR

Metals Analysis

Analyte	Result	RL	MDL	Units	DF	Prep	Analyzed	Ву	Method	Prep Method
Aluminum	21800	55	2.2	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁵
Antimony	0.39 B	2.2	0.32	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁵
Arsenic	2.0 B	2.2	0.24	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁵
Barium	113	22	0.089	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁵
Beryllium	0.21 B	0.22	0.024	-	I	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁵
Cadmium	0.14 B	0.55	0.055	mg/kg	I	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁵
Calcium	6510	550	2.1	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁵
Chromium	11.9	1.1	0.13	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁵
Cobalt	15.2	5.5	0.065	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁵
Copper	39.2	2.8	0.24	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁵
Iron	35600	110	1.8	mg/kg	2	02/01/17	02/03/17	GT	SW846 6010C ³	SW846 3050B ⁵
Lead	7.7	2.2	0.25	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁵
Magnesium	7630	550	6.5	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁵
Manganese	873	1.7	0.040	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁵
Mercury	0.031 B	0.035	0.0056	mg/kg	1	02/01/17	02/01/17	JA	SW846 7471B ¹	SW846 7471B ⁴
Nickel	5.5	4.4	0.084	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁵
Potassium	1580	1100	20	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁵
Selenium	0.51 U	2.2	0.51	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁵
Silver a	0.22 U	1.1	0.22	mg/kg	2	02/01/17	02/03/17	GT	SW846 6010C ³	SW846 3050B ⁵
Sodium	353 B	1100	4.3	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁵
Thallium	0.44 U	1.1	0.44	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁵
Vanadium	94.3	5.5	0.092	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁵
Zinc	139	5.5	0.24	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁵

(1) Instrument QC Batch: MA41277 (2) Instrument QC Batch: MA41293 (3) Instrument QC Batch: MA41305 (4) Prep QC Batch: MP98464 (5) Prep QC Batch: MP98469

(a) Elevated detection limit due to dilution required for high interfering element.





By

AC

Page 1 of 3

Client Sample ID: BKG SS-2 Lab Sample ID:

JC36372-2

Matrix:

SO - Soil

Date Sampled: 01/26/17

Date Received: 01/27/17

SW846 8270D SW846 3546

Analyzed

02/08/17

Percent Solids: 78.8

Method: Project:

BMSMC, PR

DF

1

Prep Batch

OP173

Prep Date

01/30/17

Analytical Batch E6P1601

Run #1 Run #2

Run #1

Run #2

Initial Weight

File ID

6P34797.D

Final Volume

30.5 g 1.0 ml

ABN Special List

CAS No.	Compound	Result	RL	MDL	Units	Q
95-57-8	2-Chlorophenol	ND	83	21	ug/kg	
59-50-7	4-Chloro-3-methyl phenol	ND	210	26	ug/kg	
120-83-2	2,4-Dichlorophenol	ND	210	35	ug/kg	
105-67-9	2,4-Dimethylphenol	ND	210	74	ug/kg	
51-28-5	2,4-Dinitrophenol	ND	210	160	ug/kg	
534-52-1	4,6-Dinitro-o-cresol	ND	210	45	ug/kg	
95-48-7	2-Methylphenol	ND	83	27	ug/kg	
	3&4-Methylphenol	ND	83	34	ug/kg	
88-75-5	2-Nitrophenol	ND	210	28	ug/kg	
100-02-7	4-Nitrophenol	ND	420	110	ug/kg	
87-86-5	Pentachlorophenol	ND	170	39	ug/kg	
108-95-2	Phenol	ND	83	22	ug/kg	
58-90-2	2,3,4,6-Tetrachlorophenol	ND	210	28	ug/kg	
95-95-4	2,4,5-Trichlorophenol	ND	210	31	ug/kg	
88-06-2	2,4,6-Trichlorophenol	ND	210	25	ug/kg	
83-32-9	Acenaphthene	ND	42	14	ug/kg	
208-96-8	Acenaphthylene	ND	42	21	ug/kg	
98-86-2	Acetophenone	ND	210	8.9	ug/kg	
120-12-7	Anthracene	ND	42	26	ug/kg	
1912-24-9	Atrazine	ND	83	18	ug/kg	
56-55-3	Benzo(a)anthracene	58.2	42	12	ug/kg	
205-99-2	Benzo(b)fluoranthene	77.4	42	18	ug/kg	
191-24-2	Benzo(g,h,i)perylene	38.3	42	21	ug/kg	J
207-08-9	Benzo(k)fluoranthene	ND	42	19	ug/kg	
101-55-3	4-Bromophenyl phenyl ether	ND	83	16	ug/kg	
85-68-7	Butyl benzyl phthalate	ND	83	10	ug/kg	
92-52-4	1, l'-Biphenyl	ND	83	5.7	ug/kg	
100-52-7	Benzaldehyde	52.4	210	10	ug/kg	J
91-58-7	2-Chloronaphthaiene	ND	83	9.9	ug/kg	
106-47-8	4-Chloroaniline	ND	210	15	ug/kg	
86-74-8	Carbazole	ND	83	6.0	ug/kg	
105-60-2	Caprolactam	ND	83	16	ug/kg	

fael Infante Méndez 1(= 1888

ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank



Client Sample ID: BKG SS-2 Lab Sample ID: JC36372-2

Matrix: SO - Soil Method: SW846 83

 Method:
 SW846 8270D
 SW846 3546

 Project:
 BMSMC, PR

Date Sampled: 01/26/17 **Date Received:** 01/27/17

Percent Solids: 78.8

ABN Special List

CAS No.	Compound	Result	RL	MDL	Units	Q
218-01-9	Chrysene	88.2	42	13	ug/kg	
111-91-1	bis(2-Chloroethoxy)methane	ND	83	8.9	ug/kg	
111-44-4	bis(2-Chloroethyl)ether	ND	83	18	ug/kg	
108-60-1	bis(2-Chloroisopropyl)ether	ND	83	15	ug/kg	
7005-72-3	4-Chlorophenyl phenyl ether	ND	83	13	ug/kg	
121-14-2	2,4-Dinitrotoluene	ND	42	13	ug/kg	
606-20-2	2,6-Dinitrotoluene	ND	42	21	ug/kg	
91-94-1	3,3'-Dichlorobenzidine	ND	83	35	ug/kg	
123-91-1	1,4-Dioxane	ND	42	28	ug/kg	
132-64-9	Dibenzofuran	ND	83	17	ug/kg	
84-74-2	Di-n-butyl phthalate	ND	83	6.8	ug/kg	
117-84-0	Di-n-octyl phthalate	ND	83	10	ug/kg	
84-66-2	Diethyl phthalate	ND	83	8.9	ug/kg	
131-11-3	Dimethyl phthalate	ND	83	7.4	ug/kg	4
117-81-7	bis(2-Ethylhexyl)phthalate	ND	83	9.7	ug/kg	
206-44-0	Fluoranthene	113	42	19	ug/kg	
86-73-7	Fluorene	ND	42	19	ug/kg	
118-74-1	Hexachlorobenzene	ND	83	11	ug/kg	
87-68-3	Hexachlorobutadiene	ND	42	17	ug/kg	
77-47-4	Hexachlorocyclopentadiene	ND	420	17	ug/kg	
67-72-1	Hexachloroethane	ND	210	21	ug/kg	
193-39-5	Indeno(1,2,3-cd)pyrene	37.1	42	20	ug/kg	J
78-59-1	Isophorone	ND	83	8.9	ug/kg	
90-12-0	1-Methylnaphthalene	18.7	83	8.2	ug/kg	J
91-57-6	2-Methylnaphthalene	17.8	83	9.4	ug/kg	J
88-74-4	2-Nitroaniline	ND	210	9.8	ug/kg	
99-09-2	3-Nitroaniline	ND	210	10	ug/kg	
100-01-6	4-Nitroaniline	ND	210	11	ug/kg	
91-20-3	Naphthalene	27.2	42	12	ug/kg	J
98-95-3	Nitrobenzene	ND	83	16	ug/kg	
621-64-7	N-Nitroso-di-n-propylamine	ND	83	12	ug/kg	
86-30-6	N-Nitrosodiphenylamine	ND	210	15	ug/kg	COCHOO A
85-01-8	Phenanthrene	63.1	42	14	ug/kg	age more
129-00-0	Pyrene	91.8	42	13	ug/kg	12
95-94-3	1,2,4,5-Tetrachlorobenzene	ND	210	11	ug/kg	later Infante 3
						Méndez 8
CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Lim	its	10 = 1888
367-12-4	2-Fluorophenol	55%		23-1	15%	1/100 01/100
4165-62-2	Phenol-d5	54%		27-1		CO LICENCIADO

ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

Page 3 of 3

Client Sample ID: BKG SS-2 Lab Sample ID:

JC36372-2

Date Sampled: 01/26/17

Matrix: Method: SO - Soil SW846 8270D SW846 3546 Date Received: 01/27/17 Percent Solids: 78.8

Project:

BMSMC, PR

ABN Special List

CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limits
118-79-6	2,4,6-Tribromophenol	67%		19-152%
4165-60-0	Nitrobenzene-d5	55%		26-134%
321-60-8	2-Fluorobiphenyl	59%		39-124%
1718-51-0	Terphenyl-d14	60%		36-134%



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank



Report of Analysis

By

SG

01/30/17

10-157%

Page 1 of 1

Client Sample ID: BKG SS-2 Lab Sample ID:

JC36372-2

SO - Soil

SW846 8270D BY SIM SW846 3546

Analyzed

01/30/17

Date Sampled: 01/26/17

Date Received: 01/27/17

E4M3201

Percent Solids: 78.8

OP173A

Method: Project:

Matrix:

BMSMC, PR

Prep Date **Prep Batch Analytical Batch**

Run #1 Run #2

Initial Weight

Final Volume

30.5 g

File ID

4M69658.D

Terphenyl-d14

1.0 ml

DF

1

Run #1 Run #2

1718-51-0

CAS No.	Compound	Result	RL	MDL	Units	Q
50-32-8 53-70-3	Benzo(a)pyrene Dibenzo(a,h)anthracene	39.1 7.85	4.2 4.2	1.0 0.97	ug/kg ug/kg	
CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Lim	its	
4165-60-0 321-60-8	Nitrobenzene-d5 2-Fluorobiphenyl	39% 42%			38% 48%	

65%



ND = Not detected

MDL = Method Detection Limit

J = Indicates an estimated value

RL = Reporting Limit

B = Indicates analyte found in associated method blank

E = Indicates value exceeds calibration range

Report of Analysis

Page 1 of 1

Client Sample ID: BKG SS-2 Lab Sample ID:

JC36372-2 SO - Soil

Date Sampled: Date Received:

01/26/17 01/27/17

Matrix: Method:

SW846-8015C (DAI)

Percent Solids: 78.8

Project:

BMSMC, PR

Run #1

File ID DF GH108421.D 1

Analyzed By 01/30/17 XPL **Prep Date** n/a

Prep Batch n/a

Analytical Batch GGH5640

Run #2

Run #2

CAS No.

Initial Weight

Run #1

5.2 g

Low Molecular Alcohol List

Compound Result RL MDL Units Q

64-17-5 Ethanol ND 120 84 ug/kg 78-83-1 Isobutyl Alcohol ND 120 72 ug/kg ND 120 70 ug/kg 67-63-0 Isopropyl Alcohol 120 49 71-23-8 n-Propyl Alcohol ND ug/kg 71-36-3 120 ug/kg n-Butyl Alcohol ND 66 ug/kg 78-92-2 sec-Butyl Alcohol ND 120 65 67-56-1 Methanol ND 240 58 ug/kg

CAS No. **Surrogate Recoveries** Run#1 Run# 2 Limits

52-141% 111-27-3 Hexanol 85% 111-27-3 52-141% Hexanol 87%



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

By

KD

Analyzed

02/01/17

Report of Analysis

Page 1 of 1

Client Sample ID: Lab Sample ID:

BKG SS-2 JC36372-2

SGS Accutest LabLink@945799 10:10 13-Feb-2017

File ID

4G77548.D

SO - Soil

SW846 8081B SW846 3546

Date Sampled: 01/26/17 Date Received: 01/27/17

Prep Batch

OP215

Q

Matrix: Method:

Percent Solids: 78.8

Project:

BMSMC, PR

Prep Date

01/31/17

Analytical Batch

G4G2017

Run #1 Run #2

Initial Weight

Final Volume

15.1 g

10.0 ml

DF

1

Run #1 Run #2

Pesticide TCL List

,

fael Infante Méndez

ND = Not detected

2051-24-3

2051-24-3

MDL = Method Detection Limit

128%

80%

RL = Reporting Limit

E = Indicates value exceeds calibration range

Decachlorobiphenyl

Decachlorobiphenyl

J = Indicates an estimated value

10-156%

10-156%

B = Indicates analyte found in associated method blank

Client Sample ID: BKG SS-2 Lab Sample ID: JC36372-2 Matrix: SO - Soil

Date Sampled: 01/26/17 Date Received: 01/27/17

Percent Solids: 78.8

Project:

BMSMC, PR

Metals Analysis

Analyte	Result	RL	MDL	Units	DF	Prep	Analyzed	Ву	Method	Prep Method
Aluminum	18000	63	2.5	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Antimony	0.63 B	2.5	0.37	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B 4
Arsenic	1.6 B	2.5	0.27	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B 4
Barium	98.9	25	0.10	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B 4
Beryllium	0.16 B	0.25	0.028	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B 4
Cadmium	0.16 B	0.63	0.063	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B 4
Calcium	10000	630	2.4	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Chromium	20.4	1.3	0.15	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Cobalt	13.3	6.3	0.074	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Copper	38.3	3.1	0.27	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Iron	29200	63	1.0	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Lead	9.2	2.5	0.28	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Magnesium	7230	630	7.4	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Manganese	797	1.9	0.045	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Mercury	0.050	0.033	0.0053	mg/kg	1	02/01/17	02/01/17	JA	SW846 7471B ¹	SW846 7471B ³
Nickel	13.5	5.0	0.095	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Potassium	2210	1300	23	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Selenium	0.58 U	2.5	0.58	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Silver	0.12 U	0.63	0.12	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B 4
Sodium	144 B	1300	4.9	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Thallium	0.50 U	1.3	0.50	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Vanadium	81.8	6.3	0.10	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Zinc	129	6.3	0.28	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴

(1) Instrument QC Batch: MA41277 (2) Instrument QC Batch: MA41293 (3) Prep QC Batch: MP98464 (4) Prep QC Batch: MP98469



Report of Analysis

Page 1 of 1

Client Sample ID: Lab Sample ID:

BKG SS-2 DUP

JC36372-3 SO - Soil

Matrix: Method:

SW846-8015C (DAI)

Project:

BMSMC, PR

Date Sampled: 01/26/17

Date Received: 01/27/17

Percent Solids: 77.7

Run #1

File ID DF GH108420.D 1

Analyzed 01/30/17

By **XPL** n/a

Prep Date

Prep Batch n/a

Analytical Batch GGH5640

Run #2

Initial Weight

Run #1

5.2 g

Run #2

Low Molecular Alcohol List

CAS No.	Compound	Result	RL	MDL	Units	Q
64-17-5	Ethanol	ND	120	85	ug/kg	
78-83-1 67-63-0	Isobutyl Alcohol Isopropyl Alcohol	ND ND	120 120	73 71	ug/kg ug/kg	
71-23-8 71-36-3	n-Propyl Alcohol n-Butyl Alcohol	ND ND	120 120	50 67	ug/kg ug/kg	
78-92-2 67-56-1	sec-Butyl Alcohol Methanol	ND ND	120 250	66 59	ug/kg ug/kg	
CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Lim	Limits	
111-27-3	Hexanol	101%		52 - I	41%	
111-27-3	Hexanol	104%		52-1	41%	



ND = Not detected

MDL = Method Detection Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

N = Indicates presumptive evidence of a compound

Report of Analysis

By

KD

Page 1 of 1

Client Sample ID: BKG SS-2 DUP

JC36372-3

Lab Sample ID:

Matrix:

SO - Soil

Method:

SW846 8081B SW846 3546

File ID

4G77549.D

DF

1

Date Received: 01/27/17

Date Sampled: 01/26/17

Percent Solids: 77.7

OP215

Q

Project:

BMSMC, PR

Prep Date Prep Batch

01/31/17

Analytical Batch G4G2017

Run #1 Run #2

Initial Weight

Final Volume

Analyzed

02/01/17

15.3 g

10.0 ml

Run #1 Run #2

Pesticide TCL List

CAS No.	Compound	Result	RL	MDL	Units
309-00-2	Aldrin	ND	0.84	0.40	ug/kg
319-84-6	alpha-BHC	ND	0.84	0.45	ug/kg
319-85-7	beta-BHC	ND	0.84	0.53	ug/kg
319-86-8	delta-BHC	ND	0.84	0.38	ug/kg
58-89-9	gamma-BHC (Lindane)	ND	0.84	0.37	ug/kg
5103-71-9	alpha-Chlordane	ND	0.84	0.40	ug/kg
5103-74-2	gamma-Chlordane	ND	0.84	0.37	ug/kg
60-57-1	Dieldrin	ND	0.84	0.42	ug/kg
72-54-8	4,4'-DDD	ND	0.84	0.54	ug/kg
72-55-9	4,4'-DDE	ND	0.84	0.44	ug/kg
50-29-3	4,4'-DDT	ND	0.84	0.50	ug/kg
72-20-8	Endrin	ND	0.84	0.39	ug/kg
1031-07-8	Endosulfan sulfate	ND	0.84	0.34	ug/kg
7421-93-4	Endrin aldehyde	ND	0.84	0.50	ug/kg
959-98-8	Endosulfan-I	ND	0.84	0.44	ug/kg
33213-65-9	Endosulfan-II	ND	0.84	0.44	ug/kg
76-44-8	Heptachlor	ND	0.84	0.41	ug/kg
1024-57-3	Heptachlor epoxide	ND	0.84	0.45	ug/kg
72-43-5	Methoxychlor	ND	1.7	0.42	ug/kg
53494-70-5	Endrin ketone	ND	0.84	0.65	ug/kg
8001-35-2	Toxaphene	ND	21	8.7	ug/kg
CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limi	its

Surrogate Recoveries	Run# 1	Run# 2	Limits
Tetrachloro-m-xylene	89%		24-136%
Tetrachloro-m-xylene	82%		24-136%
Decachlorobiphenyl	126%		10-153%
Decachlorobiphenyl	87%		10-153%
	Tetrachloro-m-xylene Tetrachloro-m-xylene Decachlorobiphenyl	Tetrachloro-m-xylene 89% Tetrachloro-m-xylene 82% Decachlorobiphenyl 126%	Tetrachloro-m-xylene 89% Tetrachloro-m-xylene 82% Decachlorobiphenyl 126%



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

Report of Analysis

Page 1 of 3

Client Sample ID: BKG SS-3 Lab Sample ID:

JC36372-4

SO - Soil

Date Sampled: 01/26/17 Date Received:

01/27/17

Matrix: Method:

SW846 8270D SW846 3546

Percent Solids: 86.7

Project:

BMSMC, PR

Run #1

File ID 6P34798.D DF Analyzed 02/08/17 1

Prep Date By 01/30/17 AC

Prep Batch **OP173**

Analytical Batch

Run #2

Final Volume

E6P1601

Initial Weight

30.0 g 1.0 ml

Run #1 Run #2

ABN Special List

CAS No.	Compound	Result	RL	MDL	Units	Q
95-57-8	2-Chlorophenol	ND	77	19	ug/kg	
59-50-7	4-Chloro-3-methyl phenol	ND	190	24	ug/kg	
120-83-2	2,4-Dichlorophenol	ND	190	33	ug/kg	
105-67-9	2,4-Dimethylphenol	ND	190	68	ug/kg	
51-28-5	2,4-Dinitrophenol	ND	190	140	ug/kg	
534-52-1	4,6-Dinitro-o-cresol	ND	190	41	ug/kg	
95-48-7	2-Methylphenol	ND	77	25	ug/kg	
	3&4-Methylphenol	ND	77	32	ug/kg	
88-75-5	2-Nitrophenol	ND	190	25	ug/kg	
100-02-7	4-Nitrophenol	ND	380	100	ug/kg	
87-86-5	Pentachlorophenol	ND	150	36	ug/kg	
108-95-2	Phenol	ND	77	20	ug/kg	
58-90-2	2,3,4,6-Tetrachlorophenol	ND	190	25	ug/kg	
95-95-4	2,4,5-Trichlorophenol	ND	190	29	ug/kg	
88-06-2	2,4,6-Trichlorophenol	ND	190	23	ug/kg	
83-32-9	Acenaphthene	36.6	38	13	ug/kg	J
208-96-8	Acenaphthylene	ND	38	20	ug/kg	
98-86-2	Acetophenone	18.9	190	8.3	ug/kg	J
120-12-7	Anthracene	57.8	38	24	ug/kg	
1912-24-9	Atrazine	ND	77	16	ug/kg	
56-55-3	Benzo(a)anthracene	161	38	11	ug/kg	
205-99-2	Benzo(b)fluoranthene	216	38	17	ug/kg	
191-24-2	Benzo(g,h,i)perylene	104	38	19	ug/kg	
207-08-9	Benzo(k)fluoranthene	62.4	38	18	ug/kg	
101-55-3	4-Bromophenyl phenyl ether	ND	77	15	ug/kg	
85-68-7	Butyl benzyl phthalate	ND	77	9.4	ug/kg	
92-52-4	1, I'-Biphenyl	ND	77	5.3	ug/kg	
100-52-7	Benzaldehyde	47.9	190	9.5	ug/kg	J
91-58-7	2-Chloronaphthalene	ND	77	9.2	ug/kg	
106-47-8	4-Chloroaniline	ND	190	14	ug/kg	
86-74-8	Carbazole	41.6	77	5.6	ug/kg	J
105-60-2	Caprolactam	ND	77	15	ug/kg	

tuel Infante Méndez 10 = 1888

ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

Client Sample 1D: BKG SS-3 Lab Sample ID: Matrix:

JC36372-4 SO - Soil

SW846 8270D SW846 3546

Date Sampled: 01/26/17 Date Received: 01/27/17

Method: Project:

BMSMC, PR

Percent Solids: 86.7

ABN Special List

CAS No.	Compound	Result	RL	MDL	Units	Q
218-01-9	Chrysene	168	38	12	ug/kg	
111-91-1	bis(2-Chloroethoxy)methane	ND	77	8.2	ug/kg	
111-44-4	bis(2-Chloroethyl)ether	ND	77	17	ug/kg	
108-60-1	bis(2-Chloroisopropyl)ether	ND	77	14	ug/kg	
7005-72-3	4-Chlorophenyl phenyl ether	ND	77	12	ug/kg	
121-14-2	2,4-Dinitrotoluene	ND	38	12	ug/kg	
606-20-2	2,6-Dinitrotoluene	ND	38	19	ug/kg	
91-94-1	3,3'-Dichlorobenzidine	ND	77	32	ug/kg	
123-91-1	1,4-Dioxane	ND	38	25	ug/kg	
132-64-9	Dibenzofuran	21.8	77	16	ug/kg	J
84-74-2	Di-n-butyl phthalate	ND	77	6.3	ug/kg	
117-84-0	Di-n-octyl phthalate	ND	77	9.6	ug/kg	
84-66-2	Diethyl phthalate	ND	77	8.2	ug/kg	
131-11-3	Dimethyl phthalate	ND	77	6.8	ug/kg	
117-81-7	bis(2-Ethylhexyl)phthalate	ND	77	9.0	ug/kg	
206-44-0	Fluoranthene	397	38	17	ug/kg	
86-73-7	Fluorene	22.4	38	18	ug/kg	J
118-74-1	Hexachlorobenzene	ND	77	9.7	ug/kg	
87-68-3	Hexachlorobutadiene	ND	38	15	ug/kg	
77-47-4	Hexachlorocyclopentadiene	ND	380	15	ug/kg	
67-72-1	Hexachloroethane	ND	190	19	ug/kg	
193-39-5	Indeno(1,2,3-cd)pyrene	107	38	18	ug/kg	
78-59-1	Isophorone	ND	77	8.2	ug/kg	
90-12-0	1-Methylnaphthalene	ND	77	7.5	ug/kg	
91-57-6	2-Methylnaphthalene	ND	77	8.7	ug/kg	
88-74-4	2-Nitroaniline	ND	190	9.1	ug/kg	
99-09-2	3-Nitroaniline	ND	190	9.6	ug/kg	
100-01-6	4-Nitroaniline	ND	190	10	ug/kg	
91-20-3	Naphthalene	ND	38	11	ug/kg	
98-95-3	Nitrobenzene	ND	77	15	ug/kg	
621-64-7	N-Nitroso-di-n-propylamine	ND	77	11	ug/kg	
86-30-6	N-Nitrosodiphenylamine	ND	190	14	ug/kg	
85-01-8	Phenanthrene	292	38	13	ug/kg	
129-00-0	Pyrene	319	38	12	ug/kg	
95-94-3	1,2,4,5-Tetrachlorobenzene	ND	190	9.8	ug/kg	
CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Lim	its	
367-12-4	2-Fluorophenoi	74%		23-1	15%	
4165-62-2	Phenol-d5	71%		27-1	14%	



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

Page 3 of 3

Client Sample ID: BKG SS-3 Lab Sample ID:

JC36372-4 SO - Soil

Date Sampled: 01/26/17 Date Received: 01/27/17

Matrix: Method:

SW846 8270D SW846 3546

Percent Solids: 86.7

Project:

BMSMC, PR

ABN Special List

CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limits
118-79-6	2,4,6-Tribromophenol	79%		19-152%
4165-60-0	Nitrobenzene-d5	74%		26-134%
321-60-8	2-Fluorobiphenyl	72%		39-124%
1718-51-0	Terphenyl-d14	69%		36-134%





MDL = Method Detection Limit

J = Indicates an estimated value

B = Indicates analyte found in associated method blank



RL = Reporting Limit

E = Indicates value exceeds calibration range

Page 1 of 1

Client Sample ID: BKG SS-3

Lab Sample ID: Matrix:

JC36372-4

Method:

SO - Soil

SW846 8270D BY SIM SW846 3546

Date Sampled: 01/26/17

Date Received: 01/27/17

Percent Solids: 86.7

Project:

BMSMC, PR

DF Analyzed

Final Volume

By 02/02/17 SG

Prep Date 01/30/17

Prep Batch OP173A

Q

Analytical Batch

E4M3203

Run #1 Run #2

Initial Weight Run #1

File ID

4M69690.D

30.0 g

1.0 ml

1

Run #2

CAS No.	Compound	Result	RL	MDL	Units	
50-32-8 53-70-3	Benzo(a)pyrene Dibenzo(a,h)anthracene	135 23.2	3.8 3.8	0.93 0.90	ug/kg ug/kg	
CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Lim	its	
4165-60-0 321-60-8 1718-51-0	Nitrobenzene-d5 2-Fluorobiphenyl	53% 50% 72%		15-138% 12-148% 10-157%		
1/10-21-0	Terphenyl-d14	1270		10-1	2/70	

MDL = Method Detection Limit



J = Indicates an estimated value

B = Indicates analyte found in associated method blank

N = Indicates presumptive evidence of a compound

Report of Analysis

By

XPL

Page 1 of 1

Client Sample ID: BKG SS-3

Lab Sample ID:

JC36372-4 SO - Soil

Matrix: Method:

SW846-8015C (DAI)

DF

1

Date Sampled: 01/26/17

Date Received: 01/27/17

Percent Solids: 86.7

n/a

Q

Project:

BMSMC, PR

Prep Date

n/a

Analytical Batch Prep Batch GGH5640

Run #1 Run #2

Initial Weight

GH108417.D

File ID

Run #1 5.0 g

Run #2

Low Molecular Alcohol List

CAS No.	Compound	Result	RL	MDL	Units
64-17-5	Ethanol	ND	120	80	ug/kg
78-83-1	Isobutyl Alcohol	ND	120	68	ug/kg
67-63-0	Isopropyl Alcohol	ND	120	66	ug/kg
71-23-8	n-Propyl Alcohol	ND	120	46	ug/kg
71-36-3	n-Butyl Alcohol	ND	120	63	ug/kg
78-92-2	sec-Butyl Alcohol	ND	120	61	ug/kg
67-56-1	Methanol	ND	230	55	ug/kg

Analyzed

01/30/17

CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limits
111-27-3	Hexanol	76%		52-141%
111-27-3	Hexanol	71%		52-141%



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

Page 1 of 1

Client Sample ID: BKG SS-3

Lab Sample ID: JC36372-4

Matrix: Method:

Project:

SO - Soil

SGS Accutest LabLink@945799 10:10 13-Feb-2017

SW846 8081B SW846 3546

1

Date Sampled: 01/26/17 Date Received: 01/27/17

BMSMC, PR

Percent Solids: 86.7

DF File ID

Run #1 4G77550.D Analyzed 02/01/17

By KD Prep Date 01/31/17

Prep Batch **OP215**

Q

Analytical Batch G4G2017

Run #2

Initial Weight

Final Volume

15.2 g

10.0 ml

Run #1 Run #2

Pesticide TCL List

CAS No.	Compound	Result	RL	MDL	Units
309-00-2	Aldrin	ND	0.76	0.36	ug/kg
319-84-6	alpha-BHC	ND	0.76	0.41	ug/kg
319-85-7	beta-BHC	ND	0.76	0.48	ug/kg
319-86-8	delta-BHC	ND	0.76	0.34	ug/kg
58-89-9	gamma-BHC (Lindane)	ND	0.76	0.34	ug/kg
5103-71-9	alpha-Chlordane	ND	0.76	0.36	ug/kg
5103-74-2	gamma-Chlordane	ND	0.76	0.33	ug/kg
60-57-1	Dieldrin	ND	0.76	0.38	ug/kg
72-54-8	4,4'-DDD	ND	0.76	0.49	ug/kg
72-55-9	4,4'-DDE	ND	0.76	0.40	ug/kg
50-29-3	4,4'-DDT	ND	0.76	0.45	ug/kg
72-20-8	Endrin	ND	0.76	0.36	ug/kg
1031-07-8	Endosulfan sulfate	ND	0.76	0.30	ug/kg
7421-93-4	Endrin aldehyde	ND	0.76	0.45	ug/kg
959-98-8	Endosulfan-I	ND	0.76	0.40	ug/kg
33213-65-9	Endosulfan-II	ND	0.76	0.40	ug/kg
76-44-8	Heptachlor	ND	0.76	0.37	ug/kg
1024-57-3	Heptachlor epoxide	ND	0.76	0.41	ug/kg
72-43-5	Methoxychlor	ND	1.5	0.38	ug/kg
53494-70-5	Endrin ketone	ND	0.76	0.59	ug/kg
8001-35-2	Toxaphene	ND	19	7.9	ug/kg
CASNo	Surrogata Dagovarias	Dun# 1	Dun# 2	Limite	

CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limits
877-09-8	Tetrachloro-m-xylene	84%		24-136%
877-09-8	Tetrachloro-m-xylene	78%		24-136%
2051-24-3	Decachlorobiphenyl	118%		10-153%
2051-24-3	Decachlorobiphenyl	73%		10-153%



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank



Client Sample ID: BKG SS-3 Lab Sample ID: JC36372-4 Matrix: SO - Soil

Date Sampled: 01/26/17 Date Received: 01/27/17 Percent Solids: 86.7

Project: BMSMC, PR

Metals Analysis

Analyte	Result	RL	MDL	Units	DF	Prep	Analyzed	Ву	Method	Prep Method
Aluminum	17000	55	2.2	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Antimony	0.70 B	2.2	0.32	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Arsenic	1.8 B	2.2	0.24	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Barium	97.5	22	0.090	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Beryllium	0.16 B	0.22	0.024	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Cadmium	0.17 B	0.55	0.055	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Calcium	5110	550	2.1	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Chromium	12.6	1.1	0.13	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Cobalt	10.8	5.5	0.065	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Copper	35.3	2.8	0.24	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Iron	27200	55	0.88	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Lead	29.0	2.2	0.25	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Magnesium	4820	550	6.5	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Manganese	621	1.7	0.040	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Mercury	0.10	0.037	0.0059	mg/kg	1	02/01/17	02/01/17	JA	SW846 7471B [[]	SW846 7471B ³
Nickel	6.2	4.4	0.084	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Potassium	1180	1100	20	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Selenium	0.51 U	2.2	0.51	mg/kg	1	02/01/17	02/02/17	ÐΕ	SW846 6010C ²	SW846 3050B ⁴
Silver	0.11 U	0.55	0.11	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Sodium	183 B	1100	4.3	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Thallium	0.44 U	1.1	0.44	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Vanadium	81.3	5.5	0.092	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴
Zinc	85.1	5.5	0.24	mg/kg	1	02/01/17	02/02/17	DE	SW846 6010C ²	SW846 3050B ⁴

(1) Instrument QC Batch: MA41277 (2) Instrument QC Batch: MA41293 (3) Prep QC Batch: MP98464 (4) Prep QC Batch: MP98469



Analytical Batch

E2P2932

Report of Analysis

Client Sample ID: FB012617

Lab Sample ID: JC36372-5

Matrix: Method: AQ - Field Blank Soil

SW846 8270D SW846 3510C

Date Sampled: 01/26/17 Date Received: 01/27/17

Q

Percent Solids: n/a

Project:

BMSMC, PR

SGS Accutest LabLink@945799 10:10 13-Feb-2017

File ID

DF Analyzed By **Prep Date** Prep Batch 01/31/17 OP202 2P66764.D 1 02/01/17 AC

Run #1 Run #2

> Initial Volume **Final Volume** 1.0 ml 900 ml

Run #1 Run #2

ABN Special List

CAS No.	Compound	Result	RL	MDL	Units
95-57-8	2-Chlorophenol	ND	5.6	0.91	ug/l
59-50-7	4-Chloro-3-methyl phenol	ND	5.6	0.99	ug/l
120-83-2	2,4-Dichlorophenol	ND	2.2	1.4	ug/l
105-67-9	2,4-Dimethylphenol	ND	5.6	2.7	ug/l
51-28-5	2,4-Dinitrophenol	ND	11	1.7	ug/l
534-52-1	4,6-Dinitro-o-cresol	ND	5.6	1.4	ug/l
95-48-7	2-Methylphenol	ND	2.2	0.99	ug/l
	3&4-Methylphenol	ND	2.2	0.98	ug/l
88-75-5	2-Nitrophenol	ND	5.6	1.1	ug/l
100-02-7	4-Nitrophenol	ND	11	1.3	ug/l
87-86-5	Pentachlorophenol	ND	4.4	1.5	ug/l
108-95-2	Phenol	ND	2.2	0.44	ug/l
58-90-2	2,3,4,6-Tetrachlorophenol	ND	5.6	1.6	ug/l
95-95-4	2,4,5-Trichlorophenol	ND	5.6	1.5	ug/l
88-06-2	2,4,6-Trichlorophenol	ND	5.6	1.0	ug/l
83-32-9	Acenaphthene	ND	1.1	0.21	ug/l
208-96-8	Acenaphthylene	ND	1.1	0.15	ug/l
98-86-2	Acetophenone	ND	2.2	0.23	ug/l
120-12-7	Anthracene	ND	1.1	0.23	ug/l
1912-24-9	Atrazine	ND	2.2	0.50	ug/l
100-52-7	Benzaldehyde	ND	5.6	0.32	ug/l
56-55-3	Benzo(a)anthracene	ND	1.1	0.23	ug/l
205-99-2	Benzo(b)fluoranthene	ND	1.1	0.23	ug/l
191-24-2	Benzo(g,h,i)perylene	ND	1.1	0.38	ug/l
207-08-9	Benzo(k)fluoranthene	ND	1.1	0.23	ug/l
101-55-3	4-Bromophenyl phenyl ether	ND	2.2	0.45	ug/l
85-68-7	Butyl benzyl phthalate	ND	2.2	0.51	ug/l
92-52-4	1, 1'-Biphenyl	ND	1.1	0.24	ug/l
91-58-7	2-Chloronaphthalene	ND	2.2	0.26	ug/l
106-47-8	4-Chloroaniline	ND	5.6	0.38	ug/l
86-74-8	Carbazole	ND	1.1	0.25	ug/l
105-60-2	Caprolactam	ND	2.2	0.72	ug/l

fael Infante Méndez 10 # 1888

ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank



Client Sample ID: FB012617 Lab Sample ID: JC36372-5

Matrix:

AQ - Field Blank Soil

Method: SW846 8270D SW846 3510C Date Received: 01/27/17

Date Sampled: 01/26/17

Percent Solids: n/a

Project: BMSMC, PR

ABN Special List

CAS No.	Compound	Result	RL	MDL	Units	Q
		resure	I I		Cints	¥
218-01-9	Chrysene	ND	1.1	0.20	ug/l	
111-91-1	bis(2-Chloroethoxy)methane	ND	2.2	0.31	ug/l	
111-44-4	bis(2-Chloroethyl)ether	ND	2.2	0.28	ug/l	
108-60-1	bis(2-Chloroisopropyl)ether	ND	2.2	0.45	ug/l	
7005-72-3	4-Chlorophenyl phenyl ether	ND	2.2	0.41	ug/l	
121-14-2	2,4-Dinitrotoluene	ND	1.1	0.61	ug/l	
606-20-2	2,6-Dinitrotoluene	ND	1.1	0.53	ug/l	
91-94-1	3,3'-Dichlorobenzidine	ND	2.2	0.56	ug/l	
123-91-1	1,4-Dioxane	ND	1.1	0.73	ug/l	
132-64-9	Dibenzofuran	ND	5.6	0.24	ug/l	
84-74-2	Di-n-butyl phthalate	ND	2.2	0.55	ug/l	
117-84-0	Di-n-octyl phthalate	ND	2.2	0.26	ug/l	
84-66-2	Diethyl phthalate	ND	2.2	0.29	ug/l	
131-11-3	Dimethyl phthalate	ND	2.2	0.24	ug/l	
117-81-7	bis(2-Ethylhexyl)phthalate	ND	2.2	1.8	ug/l	
206-44-0	Fluoranthene	ND	1.1	0.19	ug/l	
86-73-7	Fluorene	ND	1.1	0.19	ug/l	
118-74-1	Hexachlorobenzene	ND	1.1	0.36	ug/l	
87-68-3	Hexachlorobutadiene	ND	1.1	0.55	ug/l	
77-47-4	Hexachlorocyclopentadiene	ND	11	3.1	ug/l	
67-72-1	Hexachloroethane	ND	2.2	0.43	ug/l	
193-39-5	Indeno(1,2,3-cd)pyrene	ND	1.1	0.37	ug/l	
78-59-1	Isophorone	ND	2.2	0.31	ug/l	
90-12-0	I-Methylnaphthalene	ND	1.1	0.29	ug/l	
91-57-6	2-Methylnaphthalene	ND	1.1	0.23	ug/l	
88-74-4	2-Nitroaniline	ND	5.6	0.31	ug/l	
99-09-2	3-Nitroaniline	ND	5.6	0.43	ug/l	
100-01-6	4-Nitroaniline	ND	5.6	0.49	ug/l	
91-20-3	Naphthalene	ND	1.1	0.26	ug/l	
98-95-3	Nitrobenzene	ND	2.2	0.71	ug/l	
621-64-7	N-Nitroso-di-n-propylamine	ND	2.2	0.53	ug/l	SE ISOCIADO DE
86-30-6	N-Nitrosodiphenylamine	ND	5.6	0.25	ug/l	all all
85-01-8	Phenanthrene	ND	1.1	0.19	ug/l	
129-00-0	Pyrene	ND	1.1	0.24	ug/l	fael Infante 6
95-94-3	1,2,4,5-Tetrachlorobenzene	ND	2.2	0.41	ug/l	Méndez IC = 1888
CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Lim	its	1000
	_					MICO LICENCIADO
367-12-4	2-Fluorophenol	35%		10-1		LICENT
4165-62-2	Phenol-d5	27%		10-1	10%	



MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

Client Sample ID: FB012617

Lab Sample ID: JC36372-5

Matrix: Method:

Project:

AQ - Field Blank Soil

SW846 8270D SW846 3510C BMSMC, PR

Date Sampled: 01/26/17 Date Received: 01/27/17

Percent Solids: n/a

ABN Special List

CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limits
118-79-6	2,4,6-Tribromophenol	76%		36-151%
4165-60-0	Nitrobenzene-d5	81%		34-128%
321-60-8	2-Fluorobiphenyl	87%		38-119%
1718-51-0	Terphenyl-d14	94%		26-129%



B = Indicates analyte found in associated method blank

E = Indicates value exceeds calibration range

Report of Analysis

Page I of I

Client Sample ID: FB012617

Lab Sample ID:

JC36372-5

Matrix: Method: AQ - Field Blank Soil

DF

1

SW846 8270D BY SIM SW846 3510C

Date Received: 01/27/17

Date Sampled: 01/26/17

Percent Solids: n/a

Project:

Q

BMSMC, PR

Analyzed By 01/31/17 SG

Prep Date 01/31/17

Prep Batch OP202A

Analytical Batch E4P1143

Run #1 Run #2

Run #1

Initial Volume

Final Volume

900 ml

File ID

4P20829.D

1.0 ml

Run #2

Result

RL

0.056

Run# 2

0.11

MDL

Units

50-32-8 53-70-3

4165-60-0

321-60-8

1718-51-0

CAS No.

Benzo(a)pyrene Dibenzo(a,h)anthracene

Nitrobenzene-d5

2-Fluorobiphenyl

Terphenyl-d14

Compound

ND ND

Run# 1

0.037 ug/l 0.040 ug/l

CAS No. Surrogate Recoveries

77% 60%

72%

24-125% 19-127% 10-119%

Limits

dael Infante Méndez 1(= 1888

ND = Not detected

MDL = Method Detection Limit

J = Indicates an estimated value

RL = Reporting Limit

E = Indicates value exceeds calibration range

B = Indicates analyte found in associated method blank

Page 1 of 1

Client Sample ID: FB012617

Lab Sample ID: JC36372-5

Matrix: Method: AQ - Field Blank Soil

SW846-8015C (DAI)

Date Received: 01/27/17

Percent Solids: n/a

Date Sampled: 01/26/17

Project: BMSMC, PR

Prep Batch File ID DF Analyzed By **Prep Date XPL** Run #1 GH108409.D 01/30/17 1 n/a n/a

Analytical Batch GGH5639

Run #2

Low Molecular Alcohol List

CAS No.	Compound	Result	RL	MDL	Units	Q
64-17-5	Ethanol	ND	100	55	սջ/!	
78-83-1	Isobutyl Alcohol	ND	100	36	ug/l	
67-63-0	Isopropyl Alcohol	ND	100	68	ug/l	
71-23-8	n-Propyl Alcohol	ND	100	43	ug/l	
71-36-3	n-Butyl Alcohol	ND	100	87	ug/l	
78-92-2	sec-Butyl Alcohol	ND	100	66	ug/l	
67-56-1	Methanol	ND	200	71	ug/l	
CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Lim	Limits	
111-27-3	Hexanol	111%		56-1	45%	
111-27-3	Hexanol	94%		56-1	45%	



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

By

KD

Page 1 of 1

Client Sample ID: FB012617

Lab Sample ID:

JC36372-5

Matrix: Method: AQ - Field Blank Soil

DF

1

SW846 8081B SW846 3510C

Date Sampled: 01/26/17

Date Received: 01/27/17

Project:

Percent Solids: n/a

BMSMC, PR

Prep Date

02/01/17

Prep Batch **OP233**

Q

Analytical Batch G8G79

Run #1 Run #2

Initial Volume

File ID

8G2635.D

Final Volume

Analyzed

02/03/17

Run #1

920 ml 10.0 ml

Run #2

Pesticide TCL List

CAS No.	Compound	Result	RL	MDL	Units
309-00-2	Aldrin	ND	0.011	0.0066	ug/l
319-84-6	alpha-BHC	ND	0.011	0.0065	ug/l
319-85-7	beta-BHC	ND	0.011	0.0062	ug/l
319-86-8	delta-BHC	ND	0.011	0.0050	ug/l
58-89-9	gamma-BHC (Lindane)	ND	0.011	0.0030	ug/l
5103-71-9	alpha-Chlordane	ND	0.011	0.0050	ug/l
5103-74-2	gamma-Chlordane	ND	0.011	0.0050	ug/l
60-57-1	Dieldrin	ND	0.011	0.0039	ug/l
72-54-8	4,4'-DDD	ND	0.011	0.0041	ug/l
72-55-9	4,4'-DDE	ND	0.011	0.0067	ug/l
50-29-3	4,4'-DDT	ND	0.011	0.0054	ug/l
72-20-8	Endrin	ND	0.011	0.0055	ug/l
1031-07-8	Endosulfan sulfate	ND	0.011	0.0057	ug/l
7421-93-4	Endrin aldehyde	ND	0.011	0.0056	ug/l
53494-70-5	Endrin ketone	ND	0.011	0.0055	ug/l
959-98-8	Endosulfan-I	ND	0.011	0.0054	ug/1
33213-65-9	Endosulfan-II	ND	0.011	0.0047	ug/l
76-44-8	Heptachlor	ND	0.011	0.0041	ug/l
1024-57-3	Heptachlor epoxide	ND	0.011	0.0071	ug/l
72-43-5	Methoxychlor	ND	0.022	0.0062	ug/l
8001-35-2	Toxaphene	ND	0.27	0.20	ug/l
CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limi	ts
877-09-8	Tetrachloro-m-xylene	90%		13-13	53%
877-09-8	Tetrachloro-m-xylene	89%		13-13	53%
2051-24-3	Decachlorobiphenyl	56%		10-13	38%
2051-24-3	Decachlorobiphenyl	48%	10-13	38%	



ND = Not detected

MDL = Method Detection Limit

J = Indicates an estimated value

RL = Reporting Limit E = Indicates value exceeds calibration range B = Indicates analyte found in associated method blank



Client Sample ID: FB012617 Lab Sample ID:

JC36372-5

Matrix:

AQ - Field Blank Soil

Date Received: 01/27/17

Date Sampled: 01/26/17

Percent Solids: n/a

Project:

BMSMC, PR

Total Metals Analysis

Analyte	Result	RL	MDL	Units	DF	Prep	Analyzed	Ву	Method	Prep Method
Aluminum	21 U	200	21	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Antimony	3.3 U	6.0	3.3	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Arsenic	2.2 U	3.0	2.2	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Barium	1.6 B	200	0.44	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Beryllium	0.25 U	1.0	0.25	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Cadmium	0.40 U	3.0	0.40	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Calcium	188 B	5000	33	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Chromium	0.81 U	10	0.81	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Cobalt	0.69 U	50	0.69	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Copper	2.4 U	10	2.4	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Iron	12 U	100	12	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Lead	2.3 U	3.0	2.3	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Magnesium	85 U	5000	85	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Manganese	5.1 B	15	0.39	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Mercury	0.047 U	0.20	0.047	ug/l	1	02/01/17	02/01/17	JPM	SW846 7470A ¹	SW846 7470A ⁴
Nickel	0.76 U	10	0.76	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Potassium	120 U	10000	120	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Selenium	4.1 U	10	4.1	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Silver	0.88 U	10	0.88	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Sodium	154 B	10000	24	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Thallium	1.9 U	10	1.9	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Vanadium	0.66 U	50	0.66	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Zinc	1.3 U	20	1.3	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³

(1) Instrument QC Batch: MA41282 (2) Instrument QC Batch: MA41283 (3) Prep QC Batch: MP98461 (4) Prep QC Batch: MP98472



By

AC

01/31/17

Page 1 of 3

Client Sample ID: EB012617

Lab Sample ID: JC36372-6

Matrix: Method: AQ - Equipment Blank

DF

1

SW846 8270D SW846 3510C

Date Sampled: 01/26/17 Date Received: 01/27/17

Percent Solids: n/a

OP202

Q

E2P2932

Project:

BMSMC, PR

SGS Accutest LabLink@945799 10:10 13-Feb-2017

Analytical Batch Prep Date Prep Batch

Run #1 Run #2

> Initial Volume 900 ml

File ID

2P66765.D

Final Volume

Analyzed

02/01/17

Run #1 Run #2 1.0 ml

ABN Special List

CAS No.	Compound	Result	RL	MDL	Units
95-57-8	2-Chlorophenol	ND	5.6	0.91	ug/l
59-50-7	4-Chloro-3-methyl phenol	ND	5.6	0.99	ug/1
120-83-2	2,4-Dichlorophenol	ND	2.2	1.4	ug/l
105-67-9	2,4-Dimethylphenol	ND	5.6	2.7	ug/l
51-28-5	2,4-Dinitrophenol	ND	11	1.7	ug/l
534-52-1	4,6-Dinitro-o-cresol	ND	5.6	1.4	ug/l
95-48-7	2-Methylphenol	ND	2.2	0.99	ug/l
	3&4-Methylphenol	ND	2.2	0.98	ug/l
88-75-5	2-Nitrophenol	ND	5.6	1.1	ug/l
100-02-7	4-Nitrophenol	ND	11	1.3	ug/l
87-86-5	Pentachlorophenol	ND	4.4	1.5	ug/l
108-95-2	Phenol	ND	2.2	0.44	ug/l
58-90-2	2,3,4,6-Tetrachlorophenol	ND	5.6	1.6	ug/l
95-95-4	2,4,5-Trichlorophenol	ND	5.6	1.5	ug/l
88-06-2	2,4,6-Trichlorophenol	ND	5.6	1.0	ug/l
83-32-9	Acenaphthene	ND	1.1	0.21	ug/l
208-96-8	Acenaphthylene	ND	1.1	0.15	ug/l
98-86-2	Acetophenone	ND	2.2	0.23	ug/l
120-12-7	Anthracene	ND	1.1	0.23	ug/l
1912-24-9	Atrazine	ND	2.2	0.50	ug/l
100-52-7	Benzaldehyde	ND	5.6	0.32	ug/l
56-55-3	Benzo(a)anthracene	ND	1.1	0.23	ug/l
205-99-2	Benzo(b)fluoranthene	ND	1.1	0.23	ug/i
191-24-2	Benzo(g,h,i)perylene	ND	1.1	0.38	ug/i
207-08-9	Benzo(k)fluoranthene	ND	1.1	0.23	ug/l
101-55-3	4-Bromophenyl phenyl ether	ND	2.2	0.45	ug/l
85-68-7	Butyl benzyl phthalate	ND	2.2	0.51	սջ/1
92-52-4	1,1'-Biphenyl	ND	1.1	0.24	ug/l
91-58-7	2-Chloronaphthalene	ND	2.2	0.26	ug/l
106-47-8	4-Chloroaniline	ND	5.6	0.38	ug/l
86-74-8	Carbazole	ND	1.1	0.25	ug/l
105-60-2	Caprolactam	ND	2.2	0.72	ug/l

fael Infante 1(# 1881

ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank



Client Sample ID: EB012617 Lab Sample ID:

JC36372-6

Matrix: Method:

Project:

AQ - Equipment Blank

SW846 8270D SW846 3510C BMSMC, PR

Date Received: 01/27/17

Date Sampled: 01/26/17

Percent Solids: n/a

ABN Special List

CAS No.	Compound	Result	RL	MDL	Units	Q
218-01-9	Chrysene	ND	1.1	0.20	ug/l	
111-91-1	bis(2-Chloroethoxy)methane	ND	2.2	0.31	ug/l	
111-44-4	bis(2-Chloroethyl)ether	ND	2.2	0.28	ug/l	
108-60-1	bis(2-Chloroisopropyl)ether	ND	2.2	0.45	ug/l	
7005-72-3	4-Chlorophenyl phenyl ether	ND	2.2	0.41	ug/l	
121-14-2	2,4-Dinitrotoluene	ND	1.1	0.61	ug/l	
606-20-2	2,6-Dinitrotoluene	ND	1.1	0.53	ug/l	
91-94-1	3,3'-Dichlorobenzidine	ND	2.2	0.56	ug/l	
123-91-1	1,4-Dioxane	ND	1.1	0.73	ug/l	
132-64-9	Dibenzofuran	ND	5.6	0.24	ug/l	
84-74-2	Di-n-butyl phthalate	ND	2.2	0.55	ug/l	
117-84-0	Di-n-octyl phthalate	ND	2.2	0.26	ug/i	
84-66-2	Diethyl phthalate	ND	2.2	0.29	ug/l	
131-11-3	Dimethyl phthalate	ND	2.2	0.24	ug/l	
117-81-7	bis(2-Ethylhexyl)phthalate	ND	2.2	1.8	ug/l	
206-44-0	Fluoranthene	ND	1.1	0.19	ug/l	
86-73-7	Fluorene	ND	1.1	0.19	ug/l	
118-74-1	Hexachlorobenzene	ND	1.1	0.36	ug/l	
87-68-3	Hexachlorobutadiene	ND	1.1	0.55	ug/l	
77-47-4	Hexachlorocyclopentadiene	ND	11	3.1	ug/l	
67-72-1	Hexachloroethane	ND	2.2	0.43	ug/l	
193-39-5	Indeno(1,2,3-cd)pyrene	ND	1.1	0.37	ug/l	
78-59-1	Isophorone	ND	2.2	0.31	ug/l	
90-12-0	1-Methylnaphthalene	ND	1.1	0.29	ug/l	
91-57-6	2-Methylnaphthalene	ND	1.1	0.23	ug/l	
88-74-4	2-Nitroaniline	ND	5.6	0.31	ug/1	
99-09-2	3-Nitroaniline	ND	5.6	0.43	ug/l	
100-01-6	4-Nitroaniline	ND	5.6	0.49	ug/l	
91-20-3	Naphthalene	ND	1.1	0.26	ug/l	
98-95-3	Nitrobenzene	ND	2.2	0.71	ug/l	
621-64-7	N-Nitroso-di-n-propylamine	ND	2.2	0.53	ug/l	
86-30-6	N-Nitrosodiphenylamine	ND	5.6	0.25	ug/l	S ASOCIADO DE QUE
85-01-8	Phenanthrene	ND	1.1	0.19	ug/l	at his part of the
129-00-0	Pyrene	ND	1.1	0.24	ug/l	190
95-94-3	1,2,4,5-Tetrachlorobenzene	ND	2.2	0.41	ug/l	fael Infante
CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limi	its	Méndez
	G					(3.13)
367-12-4	2-Fluorophenol	34%		10-1	10%	The Car
4165-62-2	Phenol-d5	26%		10-1	10%	CO LICENCIPAL

ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank



Page 3 of 3

Client Sample 1D: EB012617

Lab Sample ID: JC36372-6

Matrix:

AQ - Equipment Blank

Method: Project:

SW846 8270D SW846 3510C

BMSMC, PR

Date Sampled: 01/26/17 Date Received: 01/27/17

Percent Solids: n/a

ABN Special List

CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limits
118-79-6	2,4,6-Tribromophenol	81%		36-151%
4165-60-0	Nitrobenzene-d5	88%		34-128%
321-60-8	2-Fluorobiphenyl	98%		38-119%
1718-51-0	Terphenyl-d14	100%		26-129%



B = Indicates analyte found in associated method blank

E = Indicates value exceeds calibration range

By

SG

RL

0.056

Page 1 of 1

Client Sample ID:

EB012617

Lab Sample ID:

JC36372-6

Matrix: Method: AQ - Equipment Blank

DF

1

Date Sampled: 01/26/17

Date Received: 01/27/17

Project:

Percent Solids: n/a

Q

BMSMC, PR

Analyzed

01/31/17

Prep Date 01/31/17

MDL

Prep Batch OP202A

Analytical Batch E4P1143

Run #1 Run #2

Initial Volume

Final Volume

SW846 8270D BY SIM SW846 3510C

Run #1 900 ml

File ID

4P20830.D

1.0 ml

Run #2

CAS No.	Compound
50-32-8	Benzo(a)pyrene

53-70-3 Dibenzo(a,h)anthracene

ND

Result

ND

0.11 0.040Run# 2 Limits

0.037 ug/l ug/l

24-125%

19-127%

10-119%

Units

CAS No. Surrogate Recoveries

4165-60-0 Nitrobenzene-d5 321-60-8 2-Fluorobiphenyl 1718-51-0 Terphenyl-d14

Run#1

86% 65% 76%



ND = Not detected

MDL = Method Detection Limit

J = Indicates an estimated value

RL = Reporting Limit

B = Indicates analyte found in associated method blank

E = Indicates value exceeds calibration range

Page 1 of 1

Client Sample ID: EB012617

Lab Sample ID: JC36372-6

Matrix: Method: AQ - Equipment Blank SW846-8015C (DAI)

Project: BMSMC, PR **Date Sampled:** 01/26/17

Date Received: 01/27/17 Percent Solids: n/a

	File ID	DF	Analyzed	By	Prep Date	Prep Batch	Analytical Batch
Run #1	GH108410.D	1	01/30/17	XPL	n/a	n/a	GGH5639
Run #2							

Low Molecular Alcohol List

CAS No.	Compound	Result	RL	MDL	Units	Q
64-17-5	Ethanol	ND	100	55	ug/l	
78-83-1	Isobutyl Alcohol	ND	100	36	ug/l	
67-63-0	Isopropyl Alcohol	ND	100	68	ug/l	
71-23-8	n-Propyl Alcohol	ND	100	43	ug/l	
71-36-3	n-Butyl Alcohol	ND	100	87	ug/l	
78-92-2	sec-Butyl Alcohol	ND	100	66	ug/l	
67-56-1	Methanol	ND	200	71	ug/l	
CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limits		
111-27-3	Hexanol	115%		56-145%		
111-27-3	Hexanol	102%	56-145%		45%	



MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank



SGS Accutest LabLink@945799 10:10 13-Feb-2017

Report of Analysis

Page 1 of 1

Client Sample ID: EB012617

Lab Sample ID: JC36372-6

AQ - Equipment Blank

1

Matrix: Method:

SW846 8081B SW846 3510C

Date Sampled: 01/26/17 Date Received: 01/27/17

Percent Solids: n/a

Project:

BMSMC, PR

File ID DF

Prep Date Analyzed By 02/03/17 KD 02/01/17

Prep Batch OP233

Analytical Batch G8G79

Run #1 Run #2

Initial Volume

8G2636.D

Final Volume 10.0 ml

Run #1 Run #2

910 ml

Pesticide TCL List

CAS No.	Compound	Result	RL	MDL	Units	Q
309-00-2	Aldrin	ND	0.011	0.0066	ug/l	
319-84-6	alpha-BHC	ND	0.011	0.0066	ug/l	
319-85-7	beta-BHC	ND	0.011	0.0063	ug/l	
319-86-8	delta-BHC	ND	0.011	0.0050	ug/l	
58-89-9	gamma-BHC (Lindane)	ND	0.011	0.0031	ug/l	
5103-71-9	alpha-Chlordane	ND	0.011	0.0051	ug/l	
5103-74-2	gamma-Chlordane	ND	0.011	0.0050	ug/l	
60-57-1	Dieldrin	ND	0.011	0.0040	ug/l	
72-54-8	4,4'-DDD	ND	0.011	0.0042	ug/l	
72-55-9	4,4'-DDE	ND	0.011	0.0068	ug/l	
50-29-3	4,4'-DDT	ND	0.011	0.0054	ug/l	
72-20-8	Endrin	ND	0.011	0.0055	ug/l	
1031-07-8	Endosulfan sulfate	ND	0.011	0.0058	ug/l	
7421-93-4	Endrin aldehyde	ND	0.011	0.0056	ug/l	
53494-70-5	Endrin ketone	ND	0.011	0.0056	ug/l	
959-98-8	Endosulfan-I	ND	0.011	0.0055	ug/l	
33213-65-9	Endosulfan-II	ND	0.011	0.0047	ug/l	
76-44-8	Heptachlor	ND	0.011	0.0042	ug/l	
1024-57-3	Heptachlor epoxide	ND	0.011	0.0072	ug/l	
72-43-5	Methoxychlor	ND	0.022	0.0062	ug/l	
8001-35-2	Toxaphene	ND	0.27	0.20	ug/l	
CAS No.	Surrogate Recoveries	Run# 1	Run# 2	Limi	ts	
877-09-8	Tetrachloro-m-xylene	86%		13-15	3%	
877-09-8	Tetrachloro-m-xylene	85%		13-15	3%	1
2051-24-3	Decachlorobiphenyl	45%		10-13	8%	1
2051-24-3	Decachlorobiphenyl	38%	10-138%			1



ND = Not detected

MDL = Method Detection Limit

RL = Reporting Limit

E = Indicates value exceeds calibration range

J = Indicates an estimated value

B = Indicates analyte found in associated method blank

N = Indicates presumptive evidence of a compound

Client Sample ID: EB012617

Lab Sample ID: JC36372-6

Matrix: AQ - Equipment Blank

Date Sampled: 01/26/17 **Date Received:** 01/27/17

Percent Solids: n/a

Project:

BMSMC, PR

Total Metals Analysis

Analyte	Result	RL	MDL	Units	DF	Prep	Analyzed	Ву	Method	Prep Method
Aluminum	21 U	200	21	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Antimony	3.3 U	6.0	3.3	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Arsenic	2.2 U	3.0	2.2	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Barium	1.9 B	200	0.44	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Beryllium	0.25 U	1.0	0.25	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Cadmium	0.40 U	3.0	0.40	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Calcium	86.9 B	5000	33	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Chromium	0.81 U	. 10	0.81	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Cobalt	0.69 U	50	0.69	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Copper	2.4 U	10	2.4	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Iron	12 U	100	12	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Lead	2.3 U	3.0	2.3	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Magnesium	85 U	5000	85	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Manganese	2.1 B	15	0.39	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Mercury	0.067 B	0.20	0.047	ug/l	1	02/01/17	02/01/17	JPM	SW846 7470A ¹	SW846 7470A ⁴
Nickel	1.0 B	10	0.76	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Potassium	120 U	10000	120	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Selenium	4.1 U	10	4.1	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Silver	0.88 U	10	0.88	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Sodium	65.7 B	10000	24	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Thallium	1.9 U	10	1.9	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Vanadium	0.66 U	50	0.66	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³
Zinc	1.3 U	20	1.3	ug/l	1	01/31/17	02/01/17	DE	SW846 6010C ²	SW846 3010A ³

(1) Instrument QC Batch: MA41282(2) Instrument QC Batch: MA41283(3) Prep QC Batch: MP98461(4) Prep QC Batch: MP98472



Page 1 of 3

Method: SW846 8270D

Job Number: JC36372

Account:

AMANYWP Anderson, Mulholland & Associates

Project: BMSMC, PR

Sample	File ID	DF	Analyzed	By	Prep Date	Prep Batch	Analytical Batch
OP173-MS	6P34799.D	I	02/08/17	AC	01/30/17	OP173	E6P1601
OP173-MSD	6P34800.D	1	02/08/17	AC	01/30/17	OP173	E6P1601
JC36372-4	6P34798.D	1	02/08/17	AC	01/30/17	OP173	E6P1601
1							

The QC reported here applies to the following samples:

JC36372-1, JC36372-2, JC36372-4

			2-4	Spike	MS	MS	Spike	MSD	MSD		Limits
CAS No.	Compound	ug/kg	Q	ug/kg	ug/kg	%	ug/kg	ug/kg	%	RPD	Rec/RPD
95-57-8	2-Chlorophenol	NID		1010	1200	62	1010	1120	50	7	10.127/24
59-50-7	4-Chloro-3-methyl phenol	ND ND		1910 1910	1200 1420	63 74	1910 1910	1120 1270	59 67	7 11	10-137/34 11-147/35
120-83-2	2,4-Dichlorophenol			1910	1420	76			68	11	
120-63-2		ND					1910	1300			15-140/34
	2,4-Dimethylphenol	ND		1910	1190	62	1910	1050	55	13	10-151/34
51-28-5	2,4-Dinitrophenol	ND		3820	1950	51	3820	1440	38	30	10-148/49
534-52-1	4,6-Dinitro-o-cresol	ND		1910	1090	57	1910	798	42	31	10-150/48
95-48-7	2-Methylphenol	ND		1910	1210	63	1910	1120	59	8	10-138/33
00 45 6	3&4-Methylphenol	ND		1910	1230	64	1910	1140	60	8	10-143/33
88-75-5	2-Nitrophenol	ND		1910	1540	81	1910	1460	76	5	10-150/39
100-02-7	4-Nitrophenol	ND		1910	1740	91	1910	1500	79	15	10-163/38
87-86-5	Pentachlorophenol	ND		1910	1850	97	1910	1570	82	16	10-148/39
108-95-2	Phenol	ND		1910	1240	65	1910	1150	60	8	24-114/32
58-90-2	2,3,4,6-Tetrachlorophenol	ND		1910	1430	75	1910	1250	65	13	14-140/38
95-95-4	2,4,5-Trichlorophenol	ND		1910	1470	77	1910	1290	68	13	10-146/36
88-06-2	2,4,6-Trichtorophenol	ND		1910	1460	76	1910	1310	69	11	16-148/36
83-32-9	Acenaphthene	36.6	J	1910	1550	79	1910	1370	70	12	21-136/34
208-96-8	Acenaphthylene	ND		1910	1410	74	1910	1270	67	10	10-143/36
98-86-2	Acetophenone	18.9	J	1910	1310	68	1910	1260	65	4	24-127/31
120-12-7	Anthracene	57.8		1910	1480	74	1910	1340	67	10	10-147/39
1912-24-9	Atrazine	ND		1910	1440	75	1910	1310	69	9	10-161/38
56-55-3	Benzo(a)anthracene	161		1910	1670	79	1910	1420	66	16	10-151/41
50-32-8	Benzo(a)pyrene	143		1910	1590	76	1910	1410	66	12	10-149/40
205-99-2	Benzo(b)fluoranthene	216		1910	1590	72	1910	1420	63	11	10-147/42
191-24-2	Benzo(g,h,i)perylene	104		1910	1700	84	1910	1550	76	9	10-150/41
207-08-9	Benzo(k)fluoranthene	62.4		1910	1440	72	1910	1310	65	9	12-142/41
101-55-3	4-Bromophenyl phenyl ether	ND		1910	1570	82	1910	1380	72	13	26-138/37
85-68-7	Butyl benzyl phthalate	ND		1910	1610	84	1910	1400	73	14	24-143/36
92-52-4	1,1'-Biphenyl	ND		1910	1490	78	1910	1350	71	10	18-138/32
100-52-7	Benzaldehyde	47.9	J	1910	1180	59	1910	1210	61	3	10-149/37
91-58-7	2-Chloronaphthalene	ND	•	1910	1520	80	1910	1390	73	9	24-130/31
106-47-8	4-Chloroaniline	ND		1910	602	32	1910	453	24	28	10-111/52
86-74-8	Carbazole	41.6	J	1910	1450	74	1910	1300	66	11	12-146/39
105-60-2	Caprolactam	ND	v	1910	1160	61	1910	1010	53	14	10-147/40
218-01-9	Chrysene	168		1910	1660	78	1910	1430	66	15	10-147/40
111-91-1	bis(2-Chloroethoxy)methane	ND		1910	1460	76	1910	1330	70	9	10-131/41
111-44-4	bis(2-Chloroethyl)ether	ND			1420	74				1	12-142/35
111-44-4	ors(2-Cinorocinyi)ether	שו		1910	1420	74	1910	and the	~~	1	12-142/33

^{* =} Outside of Control Limits.



Page 2 of 3

Matrix Spike/Matrix Spike Duplicate Summary

Job Number: JC36372

Account:

AMANYWP Anderson, Mulholland & Associates

Project:

BMSMC, PR

Sample	File ID	DF	Analyzed	Ву	Prep Date	Prep Batch	Analytical Batch
OP173-MS	6P34799.D	1	02/08/17	AC	01/30/17	OP173	E6P1601
OP173-MSD	6P34800.D	1	02/08/17	AC	01/30/17	OP173	E6P1601
JC36372-4	6P34798.D	I	02/08/17	AC	01/30/17	OP173	E6P1601

The QC reported here applies to the following samples:

Method: SW846 8270D

JC36372-1, JC36372-2, JC36372-4

			2-4	Spike	MS	MS	Spike	MSD	MSD		Limits
CAS No.	Compound	ug/kg	Q	ug/kg	ug/kg	%	ug/kg	ug/kg	%	RPD	Rec/RPD
108-60-1	bis(2-Chloroisopropyl)ether	ND		1910	1310	69	1910	1280	67	2	10-137/33
7005-72-3	4-Chlorophenyl phenyl ether	ND		1910	1500	79	1910	1340	70	11	21-136/35
121-14-2	2,4-Dinitrotoluene	ND		1910	1730	91	1910	1550	81	11	14-148/41
606-20-2	2,6-Dinitrotoluene	ND		1910	1670	87	1910	1490	78	11	14-152/40
91-94-1	3,3'-Dichlorobenzidine	ND		3820	678	18	3820	667	17	2	10-137/47
123-91-1	1,4-Dioxane	ND		1910	555	29	1910	751	39	30	10-110/40
132-64-9	Dibenzofuran	21.8	J	1910	1480	76	1910	1320	68	11	17-141/36
84-74-2	Di-n-butyl phthalate	ND		1910	1380	72	1910	1260	66	9	26-137/35
117-84-0	Di-n-octyl phthalate	ND		1910	1350	71	1910	1230	64	9	23-145/36
84-66-2	Diethyl phthalate	ND		1910	1470	77	1910	1300	68	12	25-133/35
13 I-11-3	Dimethyl phthalate	ND		1910	1500	79	1910	1340	70	11	21-134/36
117-81-7	bis(2-Ethylhexyl)phthalate	ND		1910	1840	96	1910	1380	72	29	26-144/39
206-44-0	Fluoranthene	397		1910	1780	72	1910	1590	62	11	10-151/44
86-73-7	Fluorene	22.4	J	1910	1480	76	1910	1300	67	13	19-133/36
118-74-1	Hexachlorobenzene	ND		1910	1500	79	1910	1360	71	10	18-142/37
87-68-3	Hexachlorobutadiene	ND		1910	1430	75	1910	1360	71	5	16-137/32
77-47-4	Hexachlorocyclopentadiene	ND		3820	486	13	3820	375	10	26	10-150/50
67-72-1	Hexachloroethane	ND		1910	986	52	1910	915	48	7	10-131/38
193-39-5	Indeno(1,2,3-cd)pyrene	107		1910	1920	95	1910	1720	84	11	10-148/41
78-59-1	Isophorone	ND		1910	1370	72	1910	1250	65	9	11-142/33
90-12-0	1-Methylnaphthalene	ND		1910	1350	71	1910	1210	63	11	10-144/35
91-57-6	2-Methylnaphthalene	ND		1910	1450	76	1910	1330	70	9	10-141/35
88-74-4	2-Nitroaniline	ND		1910	1630	85	1910	1380	72	17	14-156/38
99-09-2	3-Nitroaniline	ND		1910	1090	57	1910	897	47	19	10-144/45
100-01-6	4-Nitroaniline	ND		1910	1040	54	1910	939	49	10	10-156/44
91-20-3	Naphthalene	ND		1910	1350	71	1910	1330	70	1	10-136/36
98-95-3	Nitrobenzene	ND		1910	1400	73	1910	1330	70	5	10-142/34
621-64-7	N-Nitroso-di-n-propylamine	ND		1910	1260	66	1910	1200	63	5	10-142/31
86-30-6	N-Nitrosodiphenylamine	ND		1910	1500	79	1910	1370	72	9	10-156/37
85-01-8	Phenanthrene	292		1910	1740	76	1910	1520	64	13	11-145/45
129-00-0	Pyrene	319		1910	1880	82	1910	1610	68	15	11-155/44
95-94-3	1,2,4,5-Tetrachlorobenzene	ND		1910	1680	88	1910	1500	79	11	23-136/32



Page 3 of 3

Matrix Spike/Matrix Spike Duplicate Summary Job Number: JC36372

Account: AMANYWP Anderson, Mulholland & Associates

Project: BMSMC, PR

Sample	File ID	DF	Analyzed	Ву	Prep Date	Prep Batch	Analytical Batch
OP173-MS	6P34799.D	1	02/08/17	AC	01/30/17	OP 173	E6P1601
OP173-MSD	6P34800.D	1	02/08/17	AC	01/30/17	OP173	E6P1601
JC36372-4	6P34798.D	1	02/08/17	AC	01/30/17	OP173	E6P1601

The QC reported here applies to the following samples:

JC36372-1, JC36372-2, JC36372-4

CAS No.	Surrogate Recoveries	MS	MSD	JC36372-4	Limits
367-12-4	2-Fluorophenol	76%	75%	74%	23-115%
4165-62-2	Phenol-d5	77%	72%	71%	27-114%
118-79-6	2,4,6-Tribromophenol	91%	83%	79%	19-152%
4165-60-0	Nitrobenzene-d5	75%	73%	74%	26-134%
321-60-8	2-Fluorobiphenyl	81%	73%	72%	39-124%
1718-51-0	Terphenyl-d14	83%	72%	69%	36-134%



Method: SW846 8270D

^{* =} Outside of Control Limits.

Matrix Spike/Matrix Spike Duplicate Summary Job Number: JC36372

Account:

AMANYWP Anderson, Mulholland & Associates

Project:

BMSMC, PR

Sample	File ID	DF	Analyzed	By	Prep Date	Prep Batch	Analytical Batch
OP173A-MS	4M69677.D	1	01/31/17	SG	01/30/17	OP173A	E4M3202
OP173A-MSD	4M69691.D	1	02/02/17	SG	01/30/17	OP173A	E4M3203
JC36372-4	4M69690.D	1	02/02/17	SG	01/30/17	OP173A	E4M3203

The QC reported here applies to the following samples:

Method: SW846 8270D BY SIM

JC36372-1, JC36372-2, JC36372-4

CAS No.	Compound	JC36372-4 ug/kg Q	Spike ug/kg	MS ug/kg	MS %	Spike ug/kg	MSD ug/kg	MSD %	RPD	Limits Rec/RPD
50-32-8 53-70-3	Benzo(a)pyrene Dibenzo(a,h)anthracene	135 23.2	37.9 37.9	220 59.9	224* a 97	37.8 37.8	564 118	1134* ^a 251* ^a		18-188/44 28-169/41
CAS No.	Surrogate Recoveries	MS	MSD	JC:	36372-4	Limits				
4165-60-0	Nitrobenzene-d5	49%	57%	53%	6	15-1389	6			
321-60-8	2-Fluorobiphenyl	61%	54%	50%	ó	12-1489	6			
1718-51-0	Terphenyl-d14	84%	74%	729	6	10-1579	6			

- (a) Outside of in house control limits due to possible sample nonhomogeneity.
- (b) Analytical precision exceeds in-house control limits.



^{* =} Outside of Control Limits.

Page 1 of 1

Matrix Spike/Matrix Spike Duplicate Summary

Job Number: JC36372

Account:

AMANYWP Anderson, Mulholland & Associates

Project: BMSMC, PR

Sample	File ID	DF	Analyzed	Ву	Prep Date	Prep Batch	Analytical Batch
JC36372-4MS	GH108418.D	1	01/30/17	XPL	n/a	n/a	GGH5640
JC36372-4MSD	GH108419.D	1	01/30/17	XPL	n/a	n/a	GGH5640
JC36372-4	GH108417.D	1	01/30/17	XPL	n/a	n/a	GGH5640

The QC reported here applies to the following samples:

Method: SW846-8015C (DAI)

JC36372-1, JC36372-2, JC36372-3, JC36372-4

CAS No.	Compound	JC36372-4 ug/kg Q	Spike ug/kg	MS ug/kg	MS %	Spike ug/kg	MSD ug/kg	MSD %	RPD	Limits Rec/RPD
64-17-5 78-83-1 67-63-0 71-23-8 71-36-3 78-92-2 67-56-1	Ethanol Isobutyl Alcohol Isopropyl Alcohol n-Propyl Alcohol n-Butyl Alcohol sec-Butyl Alcohol Methanol	ND ND ND ND ND ND	5770 5770 5770 5770 5770 5770 5770	6290 6650 6580 6490 6730 7450 6430	109 115 114 113 117 129 111	5770 5770 5770 5770 5770 5770 5770	4730 5550 4730 7660 6190 6690 3890	82 96 82 133 107 116 67	28* a 18 33* a 17 8 11 49* a	64-132/22 59-141/26 69-131/23 66-135/31 50-140/30 67-131/30 58-130/29
CAS No.	Surrogate Recoveries	MS	MSD	JC3	36372-4	Limits				
111-27-3 111-27-3	Hexanol Hexanol	112% 109%	93% 94%	76% 71%		52-1419 52-1419	-			

(a) High percent recoveries and no associated positive reported in the QC batch.



^{* =} Outside of Control Limits.

Matrix Spike/Matrix Spike Duplicate Summary

Job Number: JC36372

Account:

AMANYWP Anderson, Mulholland & Associates

Project:

BMSMC, PR

Sample	File ID	DF	Analyzed	By	Prep Date	Prep Batch	Analytical Batch
OP215-MS	4G77551.D	1	02/01/17	KD	01/31/17	OP215	G4G2017
OP215-MSD	4G77552.D	1	02/01/17	KD	01/31/17	OP215	G4G2017
JC36372-4	4G77550.D	1	02/01/17	KD	01/31/17	OP215	G4G2017

The QC reported here applies to the following samples:

Method: SW846 8081B

Page 1 of 1

JC36372-1, JC36372-2, JC36372-3, JC36372-4

CAS No.	Compound	JC36372-4 ug/kg Q	Spike ug/kg	MS ug/kg	MS %	Spike ug/kg	MSD ug/kg	MSD %	RPD	Limits Rec/RPD
309-00-2	Aldrin	ND	19	46.9	247* a	19	17.0	90	94* a	23-143/44
319-84-6	alpha-BHC	ND	19	17.4	92	19	17.5	92	1	18-152/47
319-85-7	beta-BHC	ND	19	9.0	47	19	16.9	89	61* a	7-143/48
319-86-8	delta-BHC	ND	19	19.5	103	19	19.8	104	2	13-155/49
58-89-9	gamma-BHC (Lindane)	ND	19	18.5	98	19	18.1	95	2	23-138/49
5103-71-9	alpha-Chlordane	ND	19	17.0	90	19	20.1	106	17	16-149/46
5103-74-2	gamma-Chlordane	ND	19	12.3	65	19	16.4	86	29	14-152/45
60-57-1	Dieldrin	ND	19	16.9	89	19	17.7	93	5	14-154/46
72-54-8	4,4' -DDD	ND	19	17.2	91	19	16.6	88	4	18-149/51
72-55-9	4,4'-DDE	ND	19	12.9	68	19	16.4	86	24	10-154/49
50-29-3	4,4'-DDT	ND	19	17.8	94	19	16.7	88	6	10-170/50
72-20-8	Endrin	ND	19	22.6	119	19	21.2	112	6	18-173/49
1031-07-8	Endosulfan sulfate	ND	19	14.3	75	19	14.8	78	3	19-132/50
7421 - 93-4	Endrin aldehyde	ND	19	25.8	136	19	16.3	86	45	10-160/53
959-98-8	Endosulfan-I	ND	19	16.3	86	19	16.7	88	2	18-143/46
33213-65-9	Endosulfan-II	ND	19	16.4	86	19	16.4	86	0	21-132/46
76-44-8	Heptachlor	ND	19	19.9	105	19	19.7	104	1	22-146/46
1024-57-3	Heptachlor epoxide	ND	19	18.0	95	19	18.9	100	5	21-151/45
72-43-5	Methoxychlor	ND	19	20.3	107	19	15.9	84	24	11-166/50
53494-70-5	Endrin ketone	ND	19	19.6	103	19	16.2	85	19	8-179/51
8001-35-2	Toxaphene	ND		ND			ND		nc	50-150/30
CAS No.	Surrogate Recoveries	MS	MSD	JC3	36372-4	Limits				
	_	· -								
877-09-8	Tetrachloro-m-xylene	83%	88%	84%		25-135%				
877-09-8	Tetrachloro-m-xylene	76%	80%			25-135%		190	CHADO D	6.
2021 24 2								OLT		A.

(a) Outside the QC limits.

112%

78%

103%

69%

118%

73%

10-156%

10-156%

Decachlorobipheny!

Decachlorobiphenyl

2051-24-3

2051-24-3





^{* =} Outside of Control Limits.

Login Number: JC36372 Account: AMANYWP - Anderson, Mulholland & Associates Project: BMSMC, PR

QC Batch ID: MP98464 Matrix Type: SOLID

Methods: SW846 7471B Units: mg/kg

Prep Date:

02/01/17

Metal	JC3637: Origina		Spikelot HGPWS1	% Rec	QC Limits	
Mercury	0.10	0.46	0.373	96.6	80-120	

Associated samples MP98464: JC36372-1, JC36372-2, JC36372-4

Results < IDL are shown as zero for calculation purposes (*) Outside of QC limits (N) Matrix Spike Rec. outside of QC limits

(anr) Analyte not requested



Login Number: JC36372 Account: AMANYWP - Anderson, Mulholland & Associates Project: BMSMC, PR

QC Batch ID: MP98464 Matrix Type: SOLID

Methods: SW846 7471B Units: mg/kg

Prep Date:

02/01/17

Metal	JC36373 Origina		Spikelot HGPWS1	% Rec	MSD RPD	QC Limit
Mercury	0.10	0.46	0.378	95.1	0.0	20

Associated samples MP98464: JC36372-1, JC36372-2, JC36372-4

Results < IDL are shown as zero for calculation purposes (*) Outside of QC limits

(N) Matrix Spike Rec. outside of QC limits

(anr) Analyte not requested



Login Number: JC36372 Account: AMANYWP - Anderson, Mulholland & Associates Project: BMSMC, PR

QC Batch ID: MP98469 Matrix Type: SOLID

Methods: SW846 6010C Units: mg/kg

Prep Date: 02/01/17

.4.

Prep Date:				02/01/17	
Metal	JC36372 Origina		Spikelo MPSPK2	t % Rec	QC Limits
Aluminum	17000	31800	2980	513.3(a)	75-125
Antimony	0.70	103	231	44.3N(b)	75-125
Arsenic	1.8	247	231	106.3	75-125
Barium	97.5	437	231	147.2N(b	75-125
Beryllium	0.16	257	231	111.3	75-125
Bismuth					
Boron					
Cadmium	0.17	252	231	109.2	75-125
Calcium	5110	11600	2880	225.1N(b	75-125
Chromium	12.6	272	231	112.4	75-125
Cobalt	10.0	270	231	112.4	75-125
Copper	35.3	305	231	116.9	75-125
Iron	27200	39000	2880	409.2(a)	75-125
Lead	29.0	287	231	111.0	75-125
Lithium					
Magnesium	4820	9100	2880	140.4N(b	75-125
Manganese	621	1450	231	359.4N(b	75-125
Molybdenum	anr				
Nickel	6.2	265	231	112.2	75-125
Palladium					
Phosphorus					
Potassium	1180	4870	2880	128.0N(b	75-125
Selenium	0.0	247	231	107.1	75-125
Silicon					
Silver	0.0	30.1	28.8	104.4	75-125
Sodium	183	3550	2880	116.8	75-125
Strontium					
Sulfur	anr				
Thallium	0.21	250	231	108.3	75-125
Tin					
Titanium	anr				title linfante
Tungsten					IC = 1888
Vanadium	81.3	346	231	114.7	75-125
					MICO LICENCIA
				P:	age 1

Login Number: JC36372 Account: AMANYWP - Anderson, Mulholland & Associates

Project: BMSMC, PR

QC Batch ID: MP98469 Matrix Type: SOLID

Methods: SW846 6010C

Units: mg/kg

Prep Date:

02/01/17

Metal	JC36372- Original		Spikelot MPSPK2	% Rec	QC Limits	
Zinc	85.1	339	231	110.1	75-125	

Zirconium

Associated samples MP98469: JC36372-1, JC36372-2, JC36372-4

Results < IDL are shown as zero for calculation purposes (*) Outside of QC limits (N) Matrix Spike Rec. outside of QC limits

- (anr) Analyte not requested

 (a) Spike amount low relative to the sample amount. Refer to lab control or spike blank for recovery
- (b) Spike recovery indicates possible matrix interference and/or sample nonhomogeneity.



Login Number: JC36372 Account: AMANYWP - Anderson, Mulholland & Associates Project: BMSMC, PR

QC Batch ID: MP98469 Matrix Type: SOLID

Methods: SW846 6010C

Units: mg/kg

Prep Date:

02/01/17

Prep Date:					02/01/1	.7	
Metal	JC36372 Origina		Spikelo MPSPK2	t % Rec	MSD RPD		QC Limit
Aluminum	17000	25900	2800	317.9(a)	20.5 (b)	20
Antimony	0.70	103	224	45.7N(c)	0.0		20
Arsenic	1.8	228	224	101.0	8.0		20
Barium	97.5	338	224	107.4	25.5 (b)	20
Beryllium	0.16	237	224	105.8	8.1		20
Bismuth							
Boron							
Cadmium	0.17	234	224	104.4	7.4		20
Calcium	5110	0440	2800	118.9	31.5 (b)	20
Chromium	12.6	245	224	103.0	10.4		20
Cobalt	10.8	245	224	104.6	9.7		20
Copper	35.3	274	224	106.6	10.7		20
Iron	27200	29800	2800	92.9	26.7 (b)	20
Lead	29.0	256	224	101.4	11.4		20
Lithium							
Magnesium	4820	7540	2800	97.2	18.8		20
Manganese	621	848	224	101.4	52,4 (b)	20
Molybdenum	anr						
Nickel	6.2	242	224	105.3	9.1		20
Palladium							
Phosphorus							
Potassium	1180	4280	2800	110.7	12.9		20
Selenium	0.0	229	224	102.2	7.6		20
Silicon							
Silver	0.0	20.6	28	102.2	5.1		20
Sodium	103	3230	2800	108.8	9.4		20
Strontium							
Sulfur	anr						
Thallium	0.21	234	224	104.4	6.6		20 SE MOCHOO PE
Tin							
Titanium	anr						duel Infante) 多 Mendez
Tungsten							IC = 1888 /S
Vanadium	81.3	309	224	101.7	11.3		20
					W ==-	1	COLICENCIAD

Page 3

Login Number: JC36372 Account: AMANYWP - Anderson, Mulholland & Associates

Project: BMSMC, PR

QC Batch ID: MP98469 Matrix Type: SOLID

Methods: SW846 6010C Units: mg/kg

Prep Date:

02/01/17

Metal	JC36372- Original		Spikelot MPSPK2	% Rec	MSD RPD	QC Limit
Zinc	85.1	305	224	98.2	10.6	20

Zirconium

Associated samples MP98469: JC36372-1, JC36372-2, JC36372-4

Results < IDL are shown as zero for calculation purposes (*) Outside of QC limits

(N) Matrix Spike Rec. outside of QC limits

(anr) Analyte not requested

- (a) Spike amount low relative to the sample amount. Refer to lab control or spike blank for recovery information.
- (b) High rpd due to possible sample nonhomogeneity.
- .(c) Spike recovery indicates possible matrix interference and/or sample nonhomogeneity.



_\$	GS. A	CUTEST	NJ	50 70			GS Accu	plest - Da Dayton,	moty Fill (UK)	0					Pack to # 5641060				<u>}</u>	PAGE OF						
				W	1212	EL. 732-329		and and an	erra.										_			_	ر ل	ازسا	()	/-
HOUSE H	Client / Recy	orting Informatio	n treeters	CERSO DANO	PONTOUR AL	Project i	nforms	don	35 664	Vel 2	524	gange?	Andra	201 19	SER	eque	sted /	Litalyi	sla (a	ee TE	ST CC	DESI	teet)	NATION OF	16/23	Matrix Codes
				Project Name										ı	ĺ		- 1		- 1		- 1	- 1	}	- 1	- 1,	DW - Droking Water
Company h				Drietol Muses	Squi <u>bb Manufa</u>	cturing Co.	. Surfac	e Soll S	iampling	1						-	- 1		-	- 1		- 1		- 1	- 1	GW - Ground Water WW - Water
Andersor Seed Add	att ett	& Associates		Street	adelpo marie-		190 11 14 1	other Proplet	DATE:	r.berk	425-0	VALUE OF SHIP	COLUMN TO	Legh.		- 1	ŀ		-	le l	- 1	i	- 1			SW - Surface Water SO - Soi
1		enue, Suite 417		<u> </u>			Billing to Company	ntarmano	n (if stiffer	end fro	m Re	port to)		-						ES1	- 1	.			- 1	St Bludge
City	Şlarie)	Zo	Cet		PR	Company	k tellenen											2	FOR	ļ					SED-Sediment CI - D4
Purchasi			10577	Humação		FR	Street Ad	idi ett						\neg					Dibanz(a,h)anthracene	22		- 1		- 1		LIQ - Other Liquid IAP - Air
Prosect Cor																		8	#	COMMENTS	- 1	l			- 1	SQL - Other Solid V/P - Wice
Terry Phone #	aylor		Faxil	Cheril Purchase C)roer 8		Cey			84	p (mp		20		- [-	ξl	칠	8				- 1	- 1	FB-Field Dlank
914-2	1-0400						American							-	١,	۱	ğ	Benzo(a)pym	흥	2		- 1			-	EB-Equipment Black RB-Rinsa Blank
Sempleris	Name/E1		Phone #	Project Manager			ADMINUT	Ç.						1.	. }	=	<u> </u>	Ĕ	3	· SEE	9	ĺ	- 1		ı	TB-Trp Blank
L				Terry Taylor		Cellecton			- 1		-	d proces	d Burden			2	2		4	9	- I		- 1		- 1	
B041							Sampled			1	12	5 w 3	¥ 8		DB015LWA	5	AB\$270SL (SVDCe)	BSIM	BSIII	METALS	METAL		ŀ			LAB USE ONLY
Targer 4	Field ID / F	ont of Collector		MEDHIDI WWW.	Cane	Time	BA PA	Matru	P of builder	오물	19403	H250 HOME	P CNCO		5 5	2	₹	<u>a</u>		3	3				_	LAB USE ONLY
		55-1			1/26/17	0930	71	So	5		П	5	П		ΧÞ	< ·	$\times L$	\succ	X	X.	Δ					1)8
1		55-2			1/26/17	1030	11	50	5		П	5	\top		XIY	टाः	71	\overline{X}	X	X	X					CHRTY
1	816	33-7	N. I.O.			1030	11	50	2		П	2	+	11	र।;	7	\neg	\neg	-1							ERTA
3	BKG	55-2	PVP	 	126/17	1.	-		5		+	5	++		X		~	Z	$\overline{\mathbf{x}}$	V	7					=73
$\perp \angle \perp$	BKG	ss-3			1/26/17	1195	17	50	5		+	5	++	╁┼	X	7		ź١	$\langle \cdot \rangle$	\Q	'		-			735
911	BKG	85 - 3	MS		1/26/17	1145	15	50		-	+		++	11	<u>ڄ</u> اج	} -	H	51	×	()	5		_	-		
	BKG	SS- 3	MSD		1/26/17	1145	T	50	5	\sqcup		5	44-	11.	ΧĮ	¥		Ċ١		\sim	\sim	\vdash			-	1/11/1
5		212617			Y16/17	1145	171	FB	8	3	ı	4	$\bot \bot$		ΧÌ		7		×	15	X					
6	CR	112617)		1/26/17	1145	11	EB	8	3	į	4		11	\times	<u> </u>	<u> </u>	X	Κ.	X.	X		_			
		112411		1	7					П	Т	П	TT	П	Π,			1				IAST	A 12		MEN	15/09-
				 			_			1	1			11		T										
1				-			-	\vdash		Н	+-			 - -								LAS	el Vê	स्यान	Pill	570
				 			├ ─	-	-	╌	+		++		\rightarrow		\neg	-		-	_					
				200000000000000000000000000000000000000			to a section		Data	Defens	arat da	Informati	00	115	L/4776	777.07	75525		(37), 21	Com	ments.	Specii	i Instru	CROTES :	905-	TRUMA INTO ALL
7.503500	Turnarou	nd Time i Busin ins di	pyN .	Amount the IRGS	Acceptant PMJ: I Date			Commer	chai "A" (L				NYASP (Categor	y A	7	Votals	to co	onsis	l of :						
15	() Std. 10 Bush	NOSE Days							cial "B" (l)		NYASP (уВ	- 1	U,\$6.	As,Ba	,8+.C	d,Ca,	Cr,Co	Fa,Pb	Mg.Mi	LI.	Sa.A	,Na,Tl.V,Zn
	5 Tlay RUSH							FULLT1 NJ Redu	(Lavel 3 ··	4}			Blata For EDD For				Incl	-00	7.	0 SI	C+V-	ine	6,,,	and C		ie lin
] 1 Day RUSH] 2 Day RUSH											_	Other			_	A	15 8	121	051	-					
] 1 Day RUSH							NJ Das	af Known			otocal Res				-										
1 5												the "B" less		OC Sun	many		Samo	la iro	ento	nu lie su	erified	unon	raceis	at in th	e Lat	oratory
Emers	payers & Rush T	A cata available VIA	Labins		mple Custody m	unt ha deen	NJ Rec	turns = R	esum • OC	. Surre amole	nary +	Part at Ra	w dalu Lasajon,	, Includ	ling co				-cinal	, .e. vi	2.2	Ma 50 475	ACCUPATION OF	TEM.	J. Sed S	CONTRACTOR STREET
Raban	Note by Sapoling	0.00	Cute Time:			IUSE DIP GIOCUI	MONEOU I	901010 00	Cit paire o	Hattra	عاليفي	d Bys	te	def	r,		118	Delle i ii	7		Resam	od Byt	(GW	12.	Lu	-
3	TWY	//Jul)	17:	26/17/20	1 -	ed E	_>			2			12	uj			1.	1/27	117	ত্য	2 Aeren	ad Bar	7		pivo	
Railrea	and by Egic	e e	Date Time:	Record By		Asimo	puistra-	d By	•					,men 1.0	اا	J ,	4	or			+1					
3 Refere	mount by:		Date Time:		Repaired By:					Cysta	dy Soc	4.10	110	m ·	tud		Preserve	1-7	Same?	77			0.5	21	Cools 2	179,20 4
5	,-				5					16	4, 11	(4) 16d	(17)	U	ent related				ட				_	- 6		-11 11 11

JC36372: Chain of Custody Page 1 of 4

EXECUTIVE NARRATIVE

SDG No:

JC36372

Laboratory:

Accutest, New Jersey

Analysis:

SW846-8270D

Number of Samples:

7

Location:

BMSMC, Humacao, PR

SUMMARY: Seven (7) samples were analyzed for selected SVOCs following method SW846-8270D and selected PAHs were also analyzed by SW846-8270D using the selective ion monitoring (SIM) technique. The sample results were assessed according to USEPA data validation guidance documents in the following order of precedence: EPA Hazardous Waste Support Section, SOP HW-35A, July 2015 –Revision 0. Semivolatile Data Validation. The QC criteria and data validation actions listed on the data review worksheets are from the primary guidance document, unless otherwise noted.

Results are valid and can be used for decision making purposes.

Critical issues:

None

Major:

None

Minor:

None

Critical findings:

None

Major findings:

None

Minor findings:

1. Initial and continuing calibration verifications meet the method and guidance document required performance criteria except in the cases described in the Data Review Worksheet. Results for were qualified as estimated (J or UJ) in affected samples.

No closing calibration verification included in data package. No action taken, professional judgment.

QC samples were not validated.

2. MS/MSD % recovery and RPD within laboratory control limits except for the cases described in the Data Review Worksheet. Results qualified as estimated (J or UJ) in sample JC36372-4 (SIM).

COMMENTS:

Results are valid and can be used for decision making purposes.

al sufant

Reviewers Name:

Rafael Infante

Chemist License 1888

Signature:

Date:

February 18, 2017

SAMPLE ORGANIC DATA SAMPLE SUMMARY

Sample ID: JC36372-1

Sample location: BMSMC, Humacao, PR

Sampling date: 1/26/2017

Matrix: Soil

Analyte Name	Result	Units [Dilution Factor	Lab Flag	Validation	Reportable
2-Chlorophenol	77	ug/kg	1	-	U	Yes
4-Chloro-3-methyl phenol	190	ug/kg	1	-	U	Yes
2,4-Dichlorophenol	190	ug/kg	1	-	U	Yes
2,4-Dimethylphenol	190	ug/kg	1	-	U	Yes
2,4-Dinitrophenol	190	ug/kg	1	-	U	Yes
4,6-Dinitro-o-cresol	190	ug/kg	1	-	U	Yes
2-Methylphenol	77	ug/kg	1	-	U	Yes
3&4-Methylphenol	77	ug/kg	1	-	U	Yes
2-Nitrophenol	190	ug/kg	1	-	U	Yes
4-Nitrophenol	390	ug/kg	1	-	U	Yes
Pentachlorophenol	150	ug/kg	1	-	U	Yes
Phenol	77	ug/kg	1	-	U	Yes
2,3,4,6-Tetrachlorophenol	190	ug/kg	1	-	UJ	Yes
2,4,5-Trichlorophenol	190	ug/kg	1	-	UJ	Yes
2,4,6-Trichlorophenol	190	ug/kg	1	-	U	Yes
Acenaphthene	48.4	ug/kg	1	-	-	Yes
Acenaphthylene	39	ug/kg	1	-	U	Yes
Acetophenone	190	ug/kg	1	-	U	Yes
Anthracene	144	ug/kg	1	-	=	Yes
Atrazine	77	ug/kg	1	-	U	Yes
Benzo(a)anthracene	628	ug/kg	1	-	-	Yes
Benzo(a)pyrene	632	ug/kg	1	-	-	Yes
Benzo(b)fluoranthene	857	ug/kg	1	-	_	Yes
Benzo(g,h,i)perylene	557	ug/kg	1	-	-	Yes
Benzo(k)fluoranthene	278	ug/kg	1	-	_	Yes
4-Bromophenyl phenyl ether	77	ug/kg	1	-	U	Yes
Butyl benzyl phthalate	77	ug/kg	1	-	U	Yes
1,1'-Biphenyl	77	ug/kg	1	-	U	Yes
Benzaldehyde	31.6	ug/kg	1	J	J	Yes
2-Chloronaphthalene	77	ug/kg	1	-	U	Yes
4-Chloroaniline	190	ug/kg	1	-	U	Yes
Carbazole	85.9	ug/kg	1	-	=	Yes
Caprolactam	77	ug/kg	1	-	U	Yes
Chrysene	620	ug/kg	1	-	-	Yes
bis(2-Chloroethoxy)methane	77	ug/kg	1	-	U	Yes

bis(2-Chloroethyl)ether	77	ug/kg	1	-	U	Yes
bis(2-Chloroisopropyl)ether	77	ug/kg	1	-	U	Yes
4-Chlorophenyl phenyl ether	77	ug/kg	1	=	U	Yes
2,4-Dinitrotoluene	39	ug/kg	1	_	U	Yes
2,6-Dinitrotoluene	39	ug/kg	1	_	U	Yes
3,3'-Dichlorobenzidine	77	ug/kg	1	-	U	Yes
1,4-Dioxane	39	ug/kg	1	-	U	Yes
Dibenzofuran	77	ug/kg	1	J	J	Yes
Di-n-butyl phthalate	77	ug/kg	1	-	U	Yes
Di-n-octyl phthalate	77	ug/kg	1	-	U	Yes
Diethyl phthalate	77	ug/kg	1	=	U	Yes
Dimethyl phthalate	77	ug/kg	1	-	U	Yes
bis (2-Ethylhexyl) phthalate	77	ug/kg	1	-	U	Yes
Fluoranthene	1250	ug/kg	1	-	-	Yes
Fluorene	34.6	ug/kg	1	J	J	Yes
Hexachlorobenzene	77	ug/kg	1	-	U	Yes
Hexachlorobutadiene	39	ug/kg	1	-	UJ	Yes
Hexachlorocyclopentadiene	390	ug/kg	1	-	U	Yes
Hexachloroethane	190	ug/kg	1	-	UJ	Yes
Indeno(1,2,3-cd)pyrene	556	ug/kg	1	-	-	Yes
Isophorone	77	ug/kg	1	-	U	Yes
1-Methylnaphthalene	77	ug/kg	1	-	U	Yes
2-Methylnaphthalene	77	ug/kg	1	-	U	Yes
2-Nitroaniline	190	ug/kg	1	-	U	Yes
3-Nitroaniline	190	ug/kg	1	-	U	Yes
4-Nitroaniline	190	ug/kg	1	-	UJ	Yes
Naphthalene	129	ug/kg	1	-	-	Yes
Nitrobenzene	77	ug/kg	1	-	U	Yes
N-Nitroso-di-n-propylamine	77	ug/kg	1	-	U	Yes
Nitrosodiphenylamine	77	ug/kg	1	-	U	Yes
Phenanthrene	661	ug/kg	1	-	U	Yes
Pyrene	1200	ug/kg	1	-	U	Yes
1,2,4,5-Tetrachlorobenzene	190	ug/kg	1	-	U	Yes
METHOD:	8270D (SIN	√ I)				

93.2

Dibenzo(a,h)anthracene

ug/kg

1

Yes

Sample location: BMSMC, Humacao, PR

Sampling date: 1/26/2017

Matrix: Soil

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
2-Chlorophenol	83	ug/kg	1	-	U	Yes
4-Chloro-3-methyl phenol	210	ug/kg	1	-	U	Yes
2,4-Dichlorophenol	210	ug/kg	1	-	U	Yes
2,4-Dimethylphenol	210	ug/kg	1	-	U	Yes
2,4-Dinitrophenol	210	ug/kg	1	-	U	Yes
4,6-Dinitro-o-cresol	210	ug/kg	1	-	U	Yes
2-Methylphenol	83	ug/kg	1	-	U	Yes
3&4-Methylphenol	83	ug/kg	1	-	U	Yes
2-Nitrophenol	210	ug/kg	1	-	U	Yes
4-Nitrophenol	420	ug/kg	1	-	U	Yes
Pentachlorophenol	170	ug/kg	1	-	U	Yes
Phenol	83	ug/kg	1	-	U	Yes
2,3,4,6-Tetrachlorophenol	210	ug/kg	1	-	UJ	Yes
2,4,5-Trichlorophenol	210	ug/kg	1	-	UJ	Yes
2,4,6-Trichlorophenol	210	ug/kg	1	-	U	Yes
Acenaphthene	42	ug/kg	1	-	U	Yes
Acenaphthylene	42	ug/kg	1	-	U	Yes
Acetophenone	210	ug/kg	1	-	U	Yes
Anthracene	42	ug/kg	1	-	U	Yes
Atrazine	83	ug/kg	1	-	U	Yes
Benzo(a)anthracene	58.2	ug/kg	1	-	-	Yes
Benzo(b)fluoranthene	77.4	ug/kg	1	-	-	Yes
Benzo(g,h,i)perylene	38.3	ug/kg	1	J	J	Yes
Benzo(k)fluoranthene	42	ug/kg	1	-	-	Yes
4-Bromophenyl phenyl ether	83	ug/kg	1	-	U	Yes
Butyl benzyl phthalate	83	ug/kg	1	-	U	Yes
1,1'-Biphenyl	83	ug/kg	1	-	U	Yes
Benzaldehyde	52.4	ug/kg	1	J	J	Yes
2-Chloronaphthalene	83	ug/kg	1	-	U	Yes
4-Chloroaniline	210	ug/kg	1	-	U	Yes
Carbazole	83	ug/kg	1	-	U	Yes
Caprolactam	83	ug/kg	1	-	U	Yes
Chrysene	88.2	ug/kg	1	-	-	Yes
bis(2-Chloroethoxy)methane	83	ug/kg	1	-	U	Yes
bis(2-Chloroethyl)ether	83	ug/kg	1	-	U	Yes
bis (2-Chlorois opropyl) ether	83	ug/kg	1	-	U	Yes
4-Chlorophenyl phenyl ether	83	ug/kg	1	-	U	Yes

2,4-Dinitrotoluene	42	ug/kg	1	=	U	Yes
2,6-Dinitrotoluene	42	ug/kg	1	=	U	Yes
3,3'-Dichlorobenzidine	83	ug/kg	1	_	U	Yes
1,4-Dioxane	42	ug/kg	1	-	U	Yes
Dibenzofuran	83	ug/kg	1	_	U	Yes
Di-n-butyl phthalate	83	ug/kg	1	_	U	Yes
Di-n-octyl phthalate	83	ug/kg	1	_	U	Yes
Diethyl phthalate	83	ug/kg	1	_	U	Yes
Dimethyl phthalate	83	ug/kg	1	=	U	Yes
bis(2-Ethylhexyl)phthalate	83	ug/kg	1	_	U	Yes
Fluoranthene	113	ug/kg	1	_	-	Yes
Fluorene	42	ug/kg	1	-	U	Yes
Hexachlorobenzene	83	ug/kg	1	_	U	Yes
Hexachlorobutadiene	42	ug/kg	1	_	UJ	Yes
Hexachlorocyclopentadiene	420	ug/kg	1	_	U	Yes
Hexachloroethane	210	ug/kg	1	_	UJ	Yes
Indeno(1,2,3-cd)pyrene	37.1	ug/kg	1	J	J	Yes
Isophorone	83	ug/kg	1	-	U	Yes
1-Methylnaphthalene	18.7	ug/kg	1	J	J	Yes
2-Methylnaphthalene	17.8	ug/kg	1	J	J	Yes
2-Nitroaniline	210	ug/kg	1	-	U	Yes
3-Nitroaniline	210	ug/kg	1	-	U	Yes
4-Nitroaniline	210	ug/kg	1	-	UJ	Yes
Naphthalene	27.2	ug/kg	1	J	J	Yes
Nitrobenzene	83	ug/kg	1	-	U	Yes
N-Nitroso-di-n-propylamine	83	ug/kg	1	-	U	Yes
Nitrosodiphenylamine	83	ug/kg	1	-	U	Yes
Phenanthrene	63.1	ug/kg	1	-	U	Yes
Pyrene	91.8	ug/kg	1	-	U	Yes
1,2,4,5-Tetrachlorobenzene	210	ug/kg	1	-	U	Yes

METHOD: 8270D (SIM)

Benzo(a)pyrene	39.1	ug/kg	1	-	-	Yes
Dibenzo(a,h)anthracene	7.85	ug/kg	1	-	-	Yes

Sample location: BMSMC, Humacao, PR

Sampling date: 1/26/2017

Matrix: Soil

	_, _,					
Analyte Name	Result	Units [Dilution Factor	Lab Flag	Validation	Reportable
2-Chlorophenol	77	ug/kg	1	-	U	Yes
4-Chloro-3-methyl phenol	190	ug/kg	1	-	U	Yes
2,4-Dichlorophenol	190	ug/kg	1	-	U	Yes
2,4-Dimethylphenol	190	ug/kg	1	-	U	Yes
2,4-Dinitrophenol	190	ug/kg	1	-	U	Yes
4,6-Dinitro-o-cresol	190	ug/kg	1	-	U	Yes
2-Methylphenol	77	ug/kg	1	-	U	Yes
3&4-Methylphenol	77	ug/kg	1	-	U	Yes
2-Nitrophenol	190	ug/kg	1	-	U	Yes
4-Nitrophenol	390	ug/kg	1	-	U	Yes
Pentachlorophenol	150	ug/kg	1	-	U	Yes
Phenol	77	ug/kg	1	-	U	Yes
2,3,4,6-Tetrachlorophenol	190	ug/kg	1	-	UJ	Yes
2,4,5-Trichlorophenol	190	ug/kg	1	-	UJ	Yes
2,4,6-Trichlorophenol	190	ug/kg	1	-	U	Yes
Acenaphthene	36.6	ug/kg	1	J	J	Yes
Acenaphthylene	38	ug/kg	1	-	U	Yes
Acetophenone	18.9	ug/kg	1	J	J	Yes
Anthracene	57.8	ug/kg	1	-	-	Yes
Atrazine	77	ug/kg	1	-	U	Yes
Benzo(a)anthracene	161	ug/kg	1	-	-	Yes
Benzo(b)fluoranthene	216	ug/kg	1	-	-	Yes
Benzo(g,h,i)perylene	104	ug/kg	1	-	-	Yes
Benzo(k)fluoranthene	62.4	ug/kg	1	-	-	Yes
4-Bromophenyl phenyl ether	77	ug/kg	1	-	U	Yes
Butyl benzyl phthalate	77	ug/kg	1	-	U	Yes
1,1'-Biphenyl	77	ug/kg	1	-	U	Yes
Benzaldehyde	47.9	ug/kg	1	J	J	Yes
2-Chloronaphthalene	77	ug/kg	1	-	U	Yes
4-Chloroaniline	190	ug/kg	1	-	U	Yes
Carbazole	41.6	ug/kg	1	J	J	Yes
Caprolactam	77	ug/kg	1	-	U	Yes
Chrysene	168	ug/kg	1	-	-	Yes
bis(2-Chloroethoxy)methane	77	ug/kg	1	-	U	Yes
bis(2-Chloroethyl)ether	77	ug/kg	1	-	U	Yes
bis(2-Chloroisopropyl)ether	77	ug/kg	1	-	U	Yes
4-Chlorophenyl phenyl ether	77	ug/kg	1	-	U	Yes

2,4-Dinitrotoluene	39	ug/kg	1	-	U	Yes
2,6-Dinitrotoluene	39	ug/kg	1	-	U	Yes
3,3'-Dichlorobenzidine	77	ug/kg	1	-	U	Yes
1,4-Dioxane	39	ug/kg	1	-	U	Yes
Dibenzofuran	21.8	ug/kg	1	J	J	Yes
Di-n-butyl phthalate	77	ug/kg	1	-	U	Yes
Di-n-octyl phthalate	77	ug/kg	1	-	U	Yes
Diethyl phthalate	77	ug/kg	1	-	U	Yes
Dimethyl phthalate	77	ug/kg	1	-	U	Yes
bis(2-Ethylhexyl)phthalate	77	ug/kg	1	-	U	Yes
Fluoranthene	397	ug/kg	1	-	-	Yes
Fluorene	22.4	ug/kg	1	J	J	Yes
Hexachlorobenzene	77	ug/kg	1	-	U	Yes
Hexachlorobutadiene	39	ug/kg	1	-	UJ	Yes
Hexachlorocyclopentadiene	390	ug/kg	1	-	U	Yes
Hexachloroethane	190	ug/kg	1	-	UJ	Yes
Indeno(1,2,3-cd)pyrene	107	ug/kg	1	-	-	Yes
Isophorone	77	ug/kg	1	-	U	Yes
1-Methylnaphthalene	77	ug/kg	1	-	U	Yes
2-Methylnaphthalene	77	ug/kg	1	-	U	Yes
2-Nitroaniline	190	ug/kg	1	-	U	Yes
3-Nitroaniline	190	ug/kg	1	-	U	Yes
4-Nitroaniline	190	ug/kg	1	-	UJ	Yes
Naphthalene	129	ug/kg	1	-	-	Yes
Nitrobenzene	77	ug/kg	1	-	U	Yes
N-Nitroso-di-n-propylamine	77	ug/kg	1	-	U	Yes
Nitrosodiphenylamine	77	ug/kg	1	-	U	Yes
Phenanthrene	292	ug/kg	1	-	U	Yes
Pyrene	319	ug/kg	1	-	U	Yes
1,2,4,5-Tetrachlorobenzene	190	ug/kg	1	-	U	Yes
METHOD	: 8270D (SIN	√ I)				
Benzo(a)pyrene	135	ug/kg	1	_	J	Yes
Dibenzo(a,h)anthracene	23.2	ug/kg	1	-	J	Yes
` ' '		5, 5				

Sample location: BMSMC, Humacao, PR

Sampling date: 1/26/2017

Matrix: AQ - Field Blank Soil

Analyte Name	Result	Units I	Dilution Factor	Lab Flag	Validation	Reportable
2-Chlorophenol	5.6	ug/l	1	-	U	Yes
4-Chloro-3-methyl phenol	5.6	ug/l	1	_	U	Yes
2,4-Dichlorophenol	2.2	ug/l	1	_	U	Yes
2,4-Dimethylphenol	5.6	ug/l	1	_	U	Yes
2,4-Dinitrophenol	11	ug/l	1	-	U	Yes
4,6-Dinitro-o-cresol	5.6	ug/l	1	-	U	Yes
2-Methylphenol	2.2	ug/l	1	-	U	Yes
3&4-Methylphenol	2.2	ug/l	1	-	U	Yes
2-Nitrophenol	5.6	ug/l	1	-	U	Yes
4-Nitrophenol	11	ug/l	1	-	U	Yes
Pentachlorophenol	4.4	ug/l	1	-	U	Yes
Phenol	2.2	ug/l	1	-	U	Yes
2,3,4,6-Tetrachlorophenol	5.6	ug/l	1	-	U	Yes
2,4,5-Trichlorophenol	5.6	ug/l	1	-	U	Yes
2,4,6-Trichlorophenol	5.6	ug/l	1	-	U	Yes
Acenaphthene	1.1	ug/l	1	-	U	Yes
Acenaphthylene	1.1	ug/l	1	-	U	Yes
Acetophenone	2.2	ug/l	1	-	U	Yes
Anthracene	1.1	ug/l	1	-	U	Yes
Atrazine	2.2	ug/l	1	-	UJ	Yes
Benzaldehyde	5.6	ug/l	1	-	U	Yes
Benzo(a)anthracene	1.1	ug/l	1	-	U	Yes
Benzo(a)pyrene	1.1	ug/l	1	-	U	Yes
Benzo(b)fluoranthene	1.1	ug/l	1	-	U	Yes
Benzo(g,h,i)perylene	1.1	ug/l	1	-	U	Yes
Benzo(k)fluoranthene	1.1	ug/l	1	-	U	Yes
4-Bromophenyl phenyl ether	1.1	ug/l	1	-	U	Yes
Butyl benzyl phthalate	2.2	ug/l	1	-	U	Yes
1,1'-Biphenyl	1.1	ug/l	1	-	U	Yes
2-Chloronaphthalene	2.2	ug/l	1	-	U	Yes
4-Chloroaniline	5.6	ug/l	1	-	U	Yes
Carbazole	1.1	ug/l	1	-	U	Yes
Caprolactam	2.2	ug/l	1	-	U	Yes
Chrysene	1.1	ug/l	1	-	U	Yes
bis(2-Chloroethoxy)methane	2.2	ug/l	1	-	UJ	Yes
bis(2-Chloroethyl)ether	2.2	ug/l	1	-	U	Yes
bis(2-Chloroisopropyl)ether	2.2	ug/l	1	-	U	Yes

4-Chlorophenyl phenyl ether	2.2	ug/l	1	-	UJ	Yes
2,4-Dinitrotoluene	1.1	ug/l	1	-	U	Yes
2,6-Dinitrotoluene	1.1	ug/l	1	-	U	Yes
3,3'-Dichlorobenzidine	2.2	ug/l	1	-	U	Yes
1,4-Dioxane	1.1	ug/l	1	-	U	Yes
Dibenzofuran	5.6	ug/l	1	-	U	Yes
Di-n-butyl phthalate	2.2	ug/l	1	-	U	Yes
Di-n-octyl phthalate	2.2	ug/l	1	-	U	Yes
Diethyl phthalate	2.2	ug/l	1	-	U	Yes
Dimethyl phthalate	2.2	ug/l	1	-	U	Yes
bis(2-Ethylhexyl)phthalate	2.2	ug/l	1	-	UJ	Yes
Fluoranthene	1.1	ug/l	1	-	U	Yes
Fluorene	1.1	ug/l	1	-	U	Yes
Hexachlorobenzene	1.1	ug/l	1	-	U	Yes
Hexachlorobutadiene	1.1	ug/l	1	-	U	Yes
Hexachlorocyclopentadiene	11	ug/l	1	-	U	Yes
Hexachloroethane	2.2	ug/l	1	-	U	Yes
Indeno(1,2,3-cd)pyrene	1.1	ug/l	1	-	U	Yes
Isophorone	2.2	ug/l	1	-	U	Yes
1-Methylnaphthalene	1.1	ug/l	1	-	U	Yes
2-Methylnaphthalene	1.1	ug/l	1	-	U	Yes
2-Nitroaniline	5.6	ug/l	1	-	UJ	Yes
3-Nitroaniline	5.6	ug/l	1	-	U	Yes
4-Nitroaniline	5.6	ug/l	1	-	U	Yes
Naphthalene	1.1	ug/l	1	-	U	Yes
Nitrobenzene	2.2	ug/l	1	-	U	Yes
N-Nitroso-di-n-propylamine	2.2	ug/l	1	-	U	Yes
Nitrosodiphenylamine	5.6	ug/l	1	-	U	Yes
Phenanthrene	1.1	ug/l	1	-	U	Yes
Pyrene	1.1	ug/l	1	-	U	Yes
1,2,4,5-Tetrachlorobenzene	2.2	ug/l	1	-	U	Yes
METHOD: 8	8270D (SIN	4)				
Benzo(a)pyrene	0.056	ug/l	1	_	U	Yes
Dibenzo(a,h)anthracene	0.11	ug/l	1	_	U	Yes
2.120(a).11/a11till accile	0.11	ч Б / і	±		•	103

Sample location: BMSMC, Humacao, PR

Sampling date: 1/26/2017

Matrix: AQ - Equipment Blank Water

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
2-Chlorophenol	5.6	ug/l	1	-	U	Yes
4-Chloro-3-methyl phenol	5.6	ug/l	1	-	U	Yes
2,4-Dichlorophenol	2.2	ug/l	1	-	U	Yes
2,4-Dimethylphenol	5.6	ug/l	1	-	U	Yes
2,4-Dinitrophenol	11	ug/l	1	-	U	Yes
4,6-Dinitro-o-cresol	5.6	ug/l	1	-	U	Yes
2-Methylphenol	2.2	ug/l	1	-	U	Yes
3&4-Methylphenol	2.2	ug/l	1	-	U	Yes
2-Nitrophenol	5.6	ug/l	1	-	U	Yes
4-Nitrophenol	11	ug/l	1	-	U	Yes
Pentachlorophenol	4.4	ug/l	1	-	U	Yes
Phenol	2.2	ug/l	1	-	U	Yes
2,3,4,6-Tetrachlorophenol	5.6	ug/l	1	-	U	Yes
2,4,5-Trichlorophenol	5.6	ug/l	1	-	U	Yes
2,4,6-Trichlorophenol	5.6	ug/l	1	-	U	Yes
Acenaphthene	1.1	ug/l	1	-	U	Yes
Acenaphthylene	1.1	ug/l	1	-	U	Yes
Acetophenone	2.2	ug/l	1	-	U	Yes
Anthracene	1.1	ug/l	1	-	U	Yes
Atrazine	2.2	ug/l	1	=	UJ	Yes
Benzaldehyde	5.6	ug/l	1	-	U	Yes
Benzo(a)anthracene	1.1	ug/l	1	-	U	Yes
Benzo(a)pyrene	1.1	ug/l	1	=	U	Yes
Benzo(b)fluoranthene	1.1	ug/l	1	=	U	Yes
Benzo(g,h,i)perylene	1.1	ug/l	1	-	U	Yes
Benzo(k)fluoranthene	1.1	ug/l	1	-	U	Yes
4-Bromophenyl phenyl ether	1.1	ug/l	1	-	U	Yes
Butyl benzyl phthalate	2.2	ug/l	1	-	U	Yes
1,1'-Biphenyl	1.1	ug/l	1	-	U	Yes
2-Chloronaphthalene	2.2	ug/l	1	-	U	Yes
4-Chloroaniline	5.6	ug/l	1	-	U	Yes
Carbazole	1.1	ug/l	1	-	U	Yes
Caprolactam	2.2	ug/l	1	-	U	Yes
Chrysene	1.1	ug/l	1	-	U	Yes
bis(2-Chloroethoxy)methane	2.2	ug/l	1	-	UJ	Yes
bis(2-Chloroethyl)ether	2.2	ug/l	1	-	U	Yes
bis(2-Chloroisopropyl)ether	2.2	ug/l	1	-	U	Yes

4-Chlorophenyl phenyl ether	2.2	ug/l	1	-	UJ	Yes
2,4-Dinitrotoluene	1.1	ug/l	1	-	U	Yes
2,6-Dinitrotoluene	1.1	ug/l	1	-	U	Yes
3,3'-Dichlorobenzidine	2.2	ug/l	1	-	U	Yes
1,4-Dioxane	1.1	ug/l	1	-	U	Yes
Dibenzofuran	5.6	ug/l	1	-	U	Yes
Di-n-butyl phthalate	2.2	ug/l	1	-	U	Yes
Di-n-octyl phthalate	2.2	ug/l	1	-	U	Yes
Diethyl phthalate	2.2	ug/l	1	-	U	Yes
Dimethyl phthalate	2.2	ug/l	1	-	U	Yes
bis(2-Ethylhexyl)phthalate	2.2	ug/l	1	-	UJ	Yes
Fluoranthene	1.1	ug/l	1	-	U	Yes
Fluorene	1.1	ug/l	1	-	U	Yes
Hexachlorobenzene	1.1	ug/l	1	-	U	Yes
Hexachlorobutadiene	1.1	ug/l	1	-	U	Yes
Hexachlorocyclopentadiene	11	ug/l	1	-	U	Yes
Hexachloroethane	2.2	ug/l	1	-	U	Yes
Indeno(1,2,3-cd)pyrene	1.1	ug/l	1	-	U	Yes
Isophorone	2.2	ug/l	1	-	U	Yes
1-Methylnaphthalene	1.1	ug/l	1	-	U	Yes
2-Methylnaphthalene	1.1	ug/l	1	-	U	Yes
2-Nitroaniline	5.6	ug/l	1	-	UJ	Yes
3-Nitroaniline	5.6	ug/l	1	-	U	Yes
4-Nitroaniline	5.6	ug/l	1	-	U	Yes
Naphthalene	1.1	ug/l	1	-	U	Yes
Nitrobenzene	1.1	ug/l	1	-	U	Yes
N-Nitroso-di-n-propylamine	2.2	ug/l	1	-	U	Yes
Nitrosodiphenylamine	2.2	ug/l	1	-	U	Yes
Phenanthrene	1.1	ug/l	1	-	U	Yes
Pyrene	1.1	ug/l	1	-	U	Yes
1,2,4,5-Tetrachlorobenzene	2.2	ug/l	1	-	U	Yes
METHOD: 8	8270D (SIN	√ I)				
Benzo(a)pyrene	0.056	ug/kg	1	-	U	Yes
Dibenzo(a,h)anthracene	0.11	ug/kg	1	-	U	Yes

Sample location: BMSMC, Humacao, PR

Sampling date: 1/26/2017

Matrix: Soil

Analyte Name	Result	Units I	Dilution Factor	Lah Flag	Validation	Renortable
2-Chlorophenol	1200	ug/kg	1	-	-	Yes
4-Chloro-3-methyl phenol	1420	ug/kg	1	_	_	Yes
2,4-Dichlorophenol	1450	ug/kg	1	_	_	Yes
2,4-Dimethylphenol	1190	ug/kg	1	_	_	Yes
2,4-Dinitrophenol	1950	ug/kg	1	_	_	Yes
4,6-Dinitro-o-cresol	1090	ug/kg	1	_	_	Yes
2-Methylphenol	1210	ug/kg	1	_	_	Yes
3&4-Methylphenol	1230	ug/kg	1	_	_	Yes
2-Nitrophenol	1540	ug/kg	1	_	_	Yes
4-Nitrophenol	1740	ug/kg	1	_	_	Yes
Pentachlorophenol	1850	ug/kg	1	_	_	Yes
Phenol	1240	ug/kg	1	_	_	Yes
2,3,4,6-Tetrachlorophenol	1430	ug/kg	1	-	-	Yes
2,4,5-Trichlorophenol	1470	ug/kg	1	_	_	Yes
2,4,6-Trichlorophenol	1460	ug/kg	1	-	-	Yes
Acenaphthene	1550	ug/kg	1	-	-	Yes
Acenaphthylene	1410	ug/kg	1	-	-	Yes
Acetophenone	1310	ug/kg	1	-	-	Yes
Anthracene	1480	ug/kg	1	-	-	Yes
Atrazine	1440	ug/kg	1	=	-	Yes
Benzo(a)anthracene	1670	ug/kg	1	=	-	Yes
Benzo(a)pyrene	1590	ug/kg	1	=	-	Yes
Benzo(b)fluoranthene	1590	ug/kg	1	-	-	Yes
Benzo(g,h,i)perylene	1700	ug/kg	1	-	-	Yes
Benzo(k)fluoranthene	1440	ug/kg	1	-	-	Yes
4-Bromophenyl phenyl ether	1570	ug/kg	1	-	-	Yes
Butyl benzyl phthalate	1610	ug/kg	1	-	-	Yes
1,1'-Biphenyl	1490	ug/kg	1	=	-	Yes
Benzaldehyde	1180	ug/kg	1	=	-	Yes
2-Chloronaphthalene	1520	ug/kg	1	-	-	Yes
4-Chloroaniline	602	ug/kg	1	-	-	Yes
Carbazole	1450	ug/kg	1	-	-	Yes
Caprolactam	1160	ug/kg	1	-	-	Yes
Chrysene	1660	ug/kg	1	-	-	Yes
bis(2-Chloroethoxy)methane	1460	ug/kg	1	-	-	Yes
bis(2-Chloroethyl)ether	1420	ug/kg	1	-	-	Yes
bis (2-Chloro is opropyl) ether	1310	ug/kg	1	-	-	Yes

4-Chlorophenyl phenyl ether	1500	ug/kg	1	-	-	Yes
2,4-Dinitrotoluene	1730	ug/kg	1	-	-	Yes
2,6-Dinitrotoluene	1670	ug/kg	1	-	-	Yes
3,3'-Dichlorobenzidine	678	ug/kg	1	-	-	Yes
1,4-Dioxane	555	ug/kg	1	-	-	Yes
Dibenzofuran	1480	ug/kg	1	-	-	Yes
Di-n-butyl phthalate	1380	ug/kg	1	-	-	Yes
Di-n-octyl phthalate	1350	ug/kg	1	-	-	Yes
Diethyl phthalate	1470	ug/kg	1	-	-	Yes
Dimethyl phthalate	1500	ug/kg	1	-	-	Yes
bis(2-Ethylhexyl)phthalate	1840	ug/kg	1	-	-	Yes
Fluoranthene	1780	ug/kg	1	-	-	Yes
Fluorene	1480	ug/kg	1	-	-	Yes
Hexachlorobenzene	1500	ug/kg	1	-	-	Yes
Hexachlorobutadiene	1430	ug/kg	1	-	-	Yes
Hexachlorocyclopentadiene	486	ug/kg	1	-	-	Yes
Hexachloroethane	986	ug/kg	1	-	-	Yes
Indeno(1,2,3-cd)pyrene	1920	ug/kg	1	-	-	Yes
Isophorone	1370	ug/kg	1	-	-	Yes
1-Methylnaphthalene	1350	ug/kg	1	-	-	Yes
2-Methylnaphthalene	1450	ug/kg	1	-	-	Yes
2-Nitroaniline	1630	ug/kg	1	-	-	Yes
3-Nitroaniline	1090	ug/kg	1	-	-	Yes
4-Nitroaniline	1040	ug/kg	1	-	-	Yes
Naphthalene	1350	ug/kg	1	-	-	Yes
Nitrobenzene	1400	ug/kg	1	-	-	Yes
N-Nitroso-di-n-propylamine	1260	ug/kg	1	-	-	Yes
Nitrosodiphenylamine	1500	ug/kg	1	-	-	Yes
Phenanthrene	1740	ug/kg	1	-	-	Yes
Pyrene	1880	ug/kg	1	-	-	Yes
1,2,4,5-Tetrachlorobenzene	1680	ug/kg	1	-	-	Yes
METHOD:	8270D (SIN	√ 1)				
Benzo(a)pyrene	220	ug/kg	1	-	-	Yes
Dibenzo(a,h)anthracene	59.9	ug/kg	1	-	-	Yes

Sample location: BMSMC, Humacao, PR

Sampling date: 1/26/2017

Matrix: AQ - Equipment Blank Soil

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
2-Chlorophenol	1120	ug/kg	1	-	-	Yes
4-Chloro-3-methyl phenol	1270	ug/kg	1	-	-	Yes
2,4-Dichlorophenol	1300	ug/kg	1	-	-	Yes
2,4-Dimethylphenol	1050	ug/kg	1	-	-	Yes
2,4-Dinitrophenol	1440	ug/kg	1	-	-	Yes
4,6-Dinitro-o-cresol	798	ug/kg	1	-	-	Yes
2-Methylphenol	1120	ug/kg	1	-	-	Yes
3&4-Methylphenol	1140	ug/kg	1	-	-	Yes
2-Nitrophenol	1460	ug/kg	1	-	-	Yes
4-Nitrophenol	1500	ug/kg	1	-	-	Yes
Pentachlorophenol	1570	ug/kg	1	-	-	Yes
Phenol	1150	ug/kg	1	-	-	Yes
2,3,4,6-Tetrachlorophenol	1250	ug/kg	1	-	-	Yes
2,4,5-Trichlorophenol	1290	ug/kg	1	=	-	Yes
2,4,6-Trichlorophenol	1310	ug/kg	1	-	-	Yes
Acenaphthene	1370	ug/kg	1	-	-	Yes
Acenaphthylene	1270	ug/kg	1	=	-	Yes
Acetophenone	1260	ug/kg	1	-	-	Yes
Anthracene	1340	ug/kg	1	=	=	Yes
Atrazine	1310	ug/kg	1	=	-	Yes
Benzo(a)anthracene	1420	ug/kg	1	=	-	Yes
Benzo(a)pyrene	1410	ug/kg	1	-	-	Yes
Benzo(b)fluoranthene	1420	ug/kg	1	=	-	Yes
Benzo(g,h,i)perylene	1550	ug/kg	1	=	-	Yes
Benzo(k)fluoranthene	1310	ug/kg	1	-	-	Yes
4-Bromophenyl phenyl ether	1380	ug/kg	1	-	-	Yes
Butyl benzyl phthalate	1400	ug/kg	1	-	-	Yes
1,1'-Biphenyl	1350	ug/kg	1	-	-	Yes
Benzaldehyde	1210	ug/kg	1	-	-	Yes
2-Chloronaphthalene	1390	ug/kg	1	-	-	Yes
4-Chloroaniline	453	ug/kg	1	-	-	Yes
Carbazole	1300	ug/kg	1	-	-	Yes
Caprolactam	1010	ug/kg	1	-	-	Yes
Chrysene	1430	ug/kg	1	-	-	Yes
bis(2-Chloroethoxy)methane	1330	ug/kg	1	-	-	Yes
bis(2-Chloroethyl)ether	1410	ug/kg	1	-	-	Yes
bis(2-Chloroisopropyl)ether	1280	ug/kg	1	-	-	Yes

4-Chlorophenyl phenyl ether	1340	ug/kg	1	-	_	Yes
2,4-Dinitrotoluene	1550	ug/kg	1	-	-	Yes
2,6-Dinitrotoluene	1490	ug/kg	1	-	-	Yes
3,3'-Dichlorobenzidine	667	ug/kg	1	-	-	Yes
1,4-Dioxane	751	ug/kg	1	-	-	Yes
Dibenzofuran	1320	ug/kg	1	-	-	Yes
Di-n-butyl phthalate	1260	ug/kg	1	-	-	Yes
Di-n-octyl phthalate	1230	ug/kg	1	-	-	Yes
Diethyl phthalate	1300	ug/kg	1	-	-	Yes
Dimethyl phthalate	1340	ug/kg	1	-	-	Yes
bis(2-Ethylhexyl)phthalate	1380	ug/kg	1	-	-	Yes
Fluoranthene	1590	ug/kg	1	-	-	Yes
Fluorene	1300	ug/kg	1	-	-	Yes
Hexachlorobenzene	1360	ug/kg	1	-	-	Yes
Hexachlorobutadiene	1360	ug/kg	1	-	-	Yes
Hexachlorocyclopentadiene	375	ug/kg	1	-	-	Yes
Hexachloroethane	915	ug/kg	1	-	-	Yes
Indeno(1,2,3-cd)pyrene	1720	ug/kg	1	-	-	Yes
Isophorone	1250	ug/kg	1	-	-	Yes
1-Methylnaphthalene	1210	ug/kg	1	-	-	Yes
2-Methylnaphthalene	1330	ug/kg	1	-	-	Yes
2-Nitroaniline	1380	ug/kg	1	-	-	Yes
3-Nitroaniline	897	ug/kg	1	-	-	Yes
4-Nitroaniline	939	ug/kg	1	-	-	Yes
Naphthalene	1330	ug/kg	1	-	-	Yes
Nitrobenzene	1330	ug/kg	1	-	-	Yes
N-Nitroso-di-n-propylamine	1200	ug/kg	1	-	-	Yes
Nitrosodiphenylamine	1370	ug/kg	1	-	-	Yes
Phenanthrene	1520	ug/kg	1	-	-	Yes
Pyrene	1610	ug/kg	1	-	-	Yes
1,2,4,5-Tetrachlorobenzene	1500	ug/kg	1	-	-	Yes
METHOD:	8270D (SIN	⁄ 1)				
Benzo(a)pyrene	564	ug/kg	1	-	_	Yes
Dibenzo(a,h)anthracene	118	ug/kg	1	-	-	Yes

	Project Number:_JC36372
	Date:January_26,_2017
	Shipping Date:January_26,_2017 EPA Region:2
4	Li A Negion2_
REVIEW OF SEMIVOLATILE	ORGANIC PACKAGE
The following guidelines for evaluating volatile or validation actions. This document will assist the rake more informed decision and in better serving results were assessed according to USEPA da following order of precedence: EPA Hazardous V 2015 – Revision 0. Semivolatile Data Validation. The Conthe data review worksheets are from the primoted.	eviewer in using professional judgment to g the needs of the data users. The sample ta validation guidance documents in the Waste Support Section, SOP HW-35A, July IC criteria and data validation actions listed
The hardcopied (laboratory name) _Accutest reviewed and the quality control and performance daincluded:	data package received has been ta summarized. The data review for SVOCs
Lab. Project/SDG No.:JC36372 No. of Samples:7_SIM/7_SCAN	Sample matrix:Soil
Trip blank No.:	
Field blank No.:JC36372-5	
Equipment blank No.:JC36372-6Field duplicate No.:	
X Data Completeness	X Laboratory Control Spikes X Field Duplicates
X Holding Times X GC/MS Tuning	X Field Duplicates
X Internal Standard Performance	X Compound Identifications
X Blanks	X Compound Quantitation
X Surrogate RecoveriesX Matrix Spike/Matrix Spike Duplicate	X Quantitation Limits
_Overall Comments:_SVOCs_TCL_special_list_analyzed _analyzed_by_method_SW846-8270D_(SIM)	d_by_method_SW846-8270D;_Selected_PAHs
Definition of Qualifiers:	
J- Estimated results	
U- Compound not detected	
R- Rejected data UJ- Estimated nondetect // /	
Calif delant	
Reviewer: February 18, 2017	
Duto ODINGIT_ 101_2011	

DATA COMPLETENESS

MISSING INFORMATION	DATE LAB. CONTACTED	<u>DATE F</u>	RECEIVED
	<u> </u>		
4			
	<u></u>	<u></u>	
	k.		
		-	
		1000	
	Y .		200
		1	<u> </u>
	- 4 4		
<u>.</u> .			<u> </u>
			1/2

All criteria were met _X_	_
Criteria were not met	
and/or see below	

HOLDING TIMES

The objective of this parameter is to ascertain the validity of the results based on the holding time of the sample from time of collection to the time of analysis.

Complete table for all samples and note the analysis and/or preservation not within criteria

SAMPLE ID	DATE SAMPLED	DATE EXTRACTED/ANALYZED	рН	ACTION
All samples extra appropriate.	acted and anal	yzed within method recomm	ende	d holding time. Sample preservation

Cooler temperature (Criteria: 4 ± 2 °C):3.6°C	
---	--

Actions

Results will be qualified based on the criteria of the following Table:

Table 1. Holding Time Actions for Semivolatile Analyses

		ing Thire Actions for Semiv	Action	
Matrix	Preserved	Criteria	Detected Associated Compounds	Non-Detected Associated Compounds
	No	≤ 7 days (for extraction) ≤ 40 days (for analysis)	Use professi	onal judgment
:	No	> 7 days (for extraction) > 40 days (for analysis)	J	Use professional judgment
Aqueous	Yes	≤ 7 days (for extraction) ≤ 40 days (for analysis)	No qualification	
	Yes	> 7 days (for extraction) > 40 days (for analysis)	J	ບJ
	Yes/No	Grossly Exceeded	J	UJ or R
	No	≤ 14 days (for extraction) ≤ 40 days (for analysis)	Use profession	onal judgment
N A	No	> 14 days (for extraction) > 40 days (for analysis)	J	Use professional judgment
Non-Aqueous Yes		≤ 14 days (for extraction) ≤ 40 days (for analysis)	No qua	lification
	Yes	> 14 days (for extraction) > 40 days (for analysis)	J	UJ
	Yes/No	Grossly Exceeded	J	UJ or R

All criteria were met __X__ Criteria were not met see below

GC/MS TUNING

The assessment of the tuning results is to determine if the sample instrumentation is within the standard tuning QC limits

- _X__ The DFTPP performance results were reviewed and found to be within the specified criteria.
- _X__ DFTPP tuning was performed for every 12 hours of sample analysis.

If no, use professional judgment to determine whether the associated data should be accepted, qualified or rejected.

Notes: These requirements do not apply when samples are analyzed by the Selected Ion Monitoring (SIM) technique.

All mass spectrometer conditions must be identical to those used during the sample analysis. Background subtraction actions resulting in spectral distortion are unacceptable

Notes: No data should be qualified based of DFTPP failure.

The requirement to analyze the instrument performance check solution is optional when analysis of PAHs/pentachlorophenol is to be performed by the SIM technique.

List	the	samples	affected:

Actions:

- 1. If sample are analyzed without a preceding valid instrument performance check or are analyzed 12 hours after the Instrument Performance Check, qualify all data in those samples as unusable (R).
- 2. If ion abundance criteria are not met, use professional judgment to determine to what extent the data may be utilized.
- 3. State in the Data Review Narrative, decisions to use analytical data associated with DFTPP instrument performance checks not meeting the contract requirements.
- 4. Use professional judgment to determine if associated data should be qualified based on the spectrum of the mass calibration compounds.

All criteria were met	_X
Criteria were not met	
and/or see below	20.00

INITIAL CALIBRATION VERIFICATION

Compliance requirements for satisfactory instrument calibration are established to ensure that the instrument is capable of producing and maintaining acceptable quantitative data.

Date of initial calibration:_ Instrument ID numbers:	01/10-11/17_(SCAN) GCMS2P	12/14/16_(SIM) GCMS4M
Matrix/Level:	Aqueous/low	Aqueous/low
	01/10/167_(SIM)	
Instrument ID numbers:	GCMS4P	GCMS6P
Matrix/Level:	Aqueous/low	Aqueous/low

DATE	LAB ID#	FILE	CRITERIA OUT RFs, %RSD, %D, r	COMPOUND	SAMPLES AFFECTED
Initial a	and init	ial calib		ts the method and guid nance criteria.	dance validation document

Note: Instruments GCMS3M; GCMS6P (02/01/17); and GCMSZ were also employed for running QC samples for this data packages. QC samples not validated.

Actions:

Qualify the initial calibration analytes listed in Table 2 using the following criteria:

Table 3. Initial Calibration Actions for Semivolatile Analysis

Cultural	Action		
Criteria	Detect	Non-detect	
Initial Calibration not performed at specified frequency and sequence	Use professional judgment	Use professional judgment R	
Initial Calibration not performed at the specified concentrations	J	UJ	
RRF < Minimum RRF in Table 2 for target analyte	Use professional judgment J+ or R	R	
RRF ≥ Minimum RRF in Table 2 for target analyte	No qualification	No qualification	
%RSD > Maximum %RSD in Table 2 for target analyte	J	Use professional judgment	
%RSD ≤ Maximum %RSD in Table 2 for target analyte	No qualification	No qualification	

Initial Calibration

 $\begin{tabular}{ll} Table 2. RRF, \% RSD, and \% D \begin{tabular}{ll} Acceptance Criteria in Initial Calibration and CCV for Semivolatile Analysis \end{tabular}$

Analyte	Minimum RRF	Maximum %RSD	Opening Maximum %D¹	Opening Maximum %D¹
1,4-Dioxane	0.010	40.0	± 40.0	± 50.0
Benzaldehyde	0.100	40.0	± 40.0	± 50.0
Phenol	0.080	20.0	±20.0	± 25.0
Bis(2-chloroethyl)ether	0.100	20.0	± 20.0	±25.0
2-Chlorophenol	0.200	20.0	± 20.0	±25.0
2-Methylphenol	0.010	20.0	±20.0	±25.0
3-Methylphenol	0.010	20.0	±20.0	± 25.0
2,2'-Oxybis-(1-chloropropane)	0.010	20.0	±25.0	± 50.0
Acetophenone	0.060	20.0	±20.0	±25.0
4-Methylphenol	0.010	20.0	±20.0	±25.0
N-Nitroso-di-n-propylamine	0.080	20.0	±25.0	±25.0
lexachloroethane	0.100	20.0	±20.0	±25.0
Nitrobenzene	0.090	20.0	± 20.0	±25.0
Isophorone	0.100	20.0	± 20.0	±25.0
2-Nitrophenol	0.060	20.0	±20.0	±25.0
2,4-Dimethylphenol	0.050	20.0	±25.0	± 50.0
Bis(2-chloroethoxy)methane	0.080	20.0	±20.0	±25.0
2,4-Dichlorophenol	0.060	20.0	±20.0	±25.0
Naphthalene	0.200	20.0	±20.0	±25.0
4-Chloroaniline	0.010	40.0	± 40.0	± 50.0
Hexachlorobutadiene	0.040	20.0	±20.0	±25.0
Caprolactam	0.010	40.0	± 30.0	± 50.0
4-Chloro-3-methylphenol	0.040	20.0	± 20.0	±25.0
2-Methylnaphthalene	0.100	20.0	±20.0	±25.0
lexachlorocyclopentadiene	0.010	40.0	± 40.0	±50.0
2,4,6-Trichlorophenol	0.090	20.0	±20.0	±25.0
2,4,5-Trichlorophenol	0.100	20.0	±20.0	±25.0
I,1'-Biphenyl	0.200	20.0	±20.0	±25.0

malyte Minimum Maximum RRF %RSD		Opening Maximum %D ¹	Opening Maximum %D ¹	
2-Chloronaphthalene	0.300	20.0	±20.0	±25.0
2-Nitroaniline	0.060	20.0	±25.0	±25.0
Dimethylphthalate	0.300	20.0	±25.0	±25.0
2,6-Dinitrotoluene	0.080	20.0	± 20.0	± 25.0
Acenaphthylene	0.400	20.0	±20.0	± 25.0
3-Nitroaniline	0.010	20.0	±25.0	± 50.0
Acenaphthene	0.200	20.0	±20.0	±25.0
2,4-Dinitrophenol	0.010	40.0	± 50.0	± 50.0
4-Nitrophenol	0.010	40.0	± 40.0	± 50.0
Dibenzofuran	0.300	20.0	±20.0	±25.0
2,4-Dinitrotoluene	0.070	20.0	±20.0	±25.0
Diethylphthalate	0.300	20.0	±20.0	±25.0
1,2,4,5-Tetrachlorobenzene	0.100	20.0	±20.0	±25.0
4-Chlorophenyl-phenylether	0.100	20.0	± 20.0	±25.0
Fluorene	0.200	20.0	±20.0	± 25.0
4-Nitroaniline	0.010	40.0	±40.0	± 50.0
4,6-Dinitro-2-methylphenol	0.010	40.0	±30.0	± 50.0
4-Bromophenyl-phenyl ether	0.070	20.0	±20.0	±25.0
N-Nitrosodiphenylamine	0.100	20.0	±20.0	±25.0
Hexachlorobenzene	0.050	20.0	±20.0	±25.0
Atrazine	0.010	40.0	±25.0	± 50.0
Pentachlorophenol	0.010	40.0	± 40.0	± 50.0
Phenanthrene	0.200	20.0	±20.0	±25.0
Anthracene	0.200	20.0	±20.0	±25.0
Carbazole	0.050	20.0	± 20.0	±25.0
Di-n-butylphthalate	0.500	20.0	± 20.0	±25.0
Fluoranthene	0.100	20.0	± 20.0	±25.0
Pyrene	0.400	20.0	±25.0	±50.0
Butylbenzylphthalate	0.100	20.0	±25.0	± 50.0

Analyte	lyte Minimum RRF		Opening Maximum %D ¹	Opening Maximum %D ¹
3,3'-Dichlorobenzidine	0.010	40.0	± 40.0	± 50.0
Benzo(a)anthracene	0.300	20.0	±20.0	±25.0
Chrysene	0.200	20.0	±20.0	± 50.0
Bis(2-ethylhexyl) phthalate	0.200	20.0	±25.0	± 50.0
Di-n-octylphthalate	0.010	40.0	± 40.0	± 50.0
Benzo(b)fluoranthene	0.010	20.0	±25.0	± 50.0
Benzo(k)fluoranthene	0.010	20.0	±25.0	± 50.0
Benzo(a)pyrene	0.010	20.0	±20.0	± 50.0
Indeno(1,2,3-cd)pyrene	0.010	20.0	±25.0	± 50.0
Dibenzo(a,h)anthracene	0.010	20.0	±25.0	± 50.0
Benzo(g,h,i)perylene	0.010	20.0	±30.0	± 50.0
2,3,4,6-Tetrachlorophenol	0.040	20.0	±20.0	± 50.0
Naphthalene	0.600	20.0	±25.0	±25.0
2-Methylnaphthalene	0.300	20.0	±20.0	±25.0
Acenaphthylene	0.900	20.0	±20.0	±25.0
Acenaphthene	0.500	20.0	±20.0	±25.0
Fluorene	0.700	20.0	±25.0	±50.0
Phenanthrene	0.300	20.0	±25.0	± 50.0
Anthracene	0.400	20.0	±25.0	± 50.0
Fluoranthene	0.400	20.0	±25.0	± 50.0
Pyrene	0.500	20.0	±30.0	± 50.0
Benzo(a)anthracene	0.400	20.0	±25.0	± 50.0
Chyrsene	0.400	20.0	±25.0	± 50.0
Benzo(b)fluoranthene	0.100	20.0	±30.0	± 50.0
Benzo(k)fluoranthene	0.100	20.0	±30.0	± 50.0
Benzo(a)pyrene	0.100	20.0	±25.0	± 50.0
Indeno(1,2,3-cd)pyrene	0.100	20.0	± 40.0	± 50.0
Dibenzo(a,h)anthracene	0.010	25.0	±40.0	± 50.0
Benzo(g,h,i)perylene	0.020	25.0	±40.0	± 50.0

Pentachlorophenol	0.010	40.0	± 50.0	±50.0		
Deuterated Monitoring Compounds						

Analyte	Minimum RRF	Maximum %RSD	Opening Maximum %D ^t	Closing Maximum %D	
1,4-Dioxane-d ₈	0.010	20.0	±25.0	±50.0	
Phenol-d ₅	0.010	20.0	±25.0	±25.0	
Bis-(2-chloroethyl)ether-ds	0.100	20.0	±20.0	±25.0	
2-Chlorophenol-d₁	0.200	20.0	±20.0	±25.0	
4-Methylphenol-d ₈	0.010	20.0	±20.0	±25.0	
4-Chloroaniline-d ₄	0.010	40.0	±40.0	±50.0	
Nitrobenzene-d ₅	0.050	20.0	±20.0	±25.0	
2-Nitrophenol-d4	0.050	20.0	±20.0	±25.0	
2,4-Dichlorophenol-d3	0.060	20.0	±20.0	±25.0	
Dimethylphthalate-d ₆	0.300	20.0	±20.0	±25.0	
Acenaphthylene-d ₈	0.400	20.0	±20.0	±25.0	
4-Nitrophenol-d4	0.010	40.0	±40.0	± 50.0	
Fluorene-d ₁₀	0.100	20.0	±20.0	±25.0	
4,6-Dinitro-2-methylphenol-d2	0.010	40.0	±30.0	± 50.0	
Anthracene-d ₁₀	0.300	20.0	± 20.0	±25.0	
Pyrene-d ₁₀	0.300	20.0	±25.0	± 50.0	
Benzo(a)pyrene-d ₁₂	0.010	20.0	± 20.0	± 50.0	
Fluoranthene-d ₁₀ (SIM)	0.400	20.0	±25.0	±50.0	
2-Methylnaphthalene-d ₁₀ (SIM)	0.300	20.0	±20.0	±25.0	

¹ If a closing CCV is acting as an opening CCV, all target analytes must meet the requirements for an opening CCV.

Note: If analysis by SIM technique is requested for PAH/pentachlorophenols, calibration standards analyzed at 0.10, 0.20, 0.40, 0.80, and 1.0 ng/uL for each target compound of interest and the associated DMCs. Pentachlorophenol will require only a four point initial calibration at 0.20, 0.40, 0.80, and 1.0 ng/uL.

Date of initial calibration:

All criteria were met _	
Criteria were not met	
and/or see below	_X

12/14/16 (SIM)

CONTINUING CALIBRATION VERIFICATION

Compliance requirements for satisfactory instrument calibration are established to ensure that the instrument is capable of producing and maintaining acceptable quantitative data.

01/10-11/17 (SCAN)

Date of co	ontinuing calibra	tion verification (CC\	/11/171 /):02/01/171	01/30/17: 02/02/17
Instrumen	it ID numbers:		GCMS2P	GCMS4M
Matrix/Lev	/el:	Aqı	ueous/low	Aqueous/low
Date of in	itial calibration v Intinuing calibra	erification (ICV):(tion verification (CCV	//17_(SCAN)	01/11/17 01/31/17: 02/01/17
Date of cle	osing CCV:	<u> </u>	6P	<u> </u>
Instrumen	t ID numbers:	GCMS6	SP	GCMS4P
Matrix/Lev	/el:	Aqı	ueous/low	Aqueous/low
DATE	LAB FILE ID#	CRITERIA OUT RFs, %RSD, <u>%D</u> , r	COMPOUND	SAMPLES AFFECTED
GCMS2P				
02/01/17	cc2898-25	30.6	1,4-dioxane*	JC36372-5; -6
		24.3	n-Nitroso-di-n-propylamine*]
		22.7	bis(2-chloroethoxy)methane	
		-44.8	2-nitroaniline]
		-21.1	2,4-dinitrophenol*]
		-28.1	4-chlorophenyl-phenyl ether]
		-26.2	4,6-dinitro-methylphenol*]
		24.5	Butylbenzylphthalate*]
		29.6	bis(2-ethylhexyl)phthalate	
02/01/17	cc2899-25	-42.0	Atrazine	
GCMS6P				
02/01/17	cc1572-50	-26.8	Hexachloroethane	JC36372-1; -2; -4
		28.0	Caprolactam*	
		-41.4	Hexachlorobutadiene]
		-20 1	2 4 5-trichlorophenol	

Note: Initial and continuing calibration verifications meet the method and guidance document required performance criteria except for the cases described in this document. Results qualified as estimated (J or UJ) in affected samples.

2,3,4,6-tetrachlorophenol 4-nitroaniline

-23.2

20.2

* % difference outside was method performance criteria but within the guidance document performance criteria. No action taken.

No action taken for QC samples.

No closing calibration verification included in data package. No action taken, professional judgment.

Actions:

Notes: Verify that the CCV is run at the required frequency (an opening and closing CCV must be run within 12-hour period).

All DMCs must meet the RRF values given in Table 2. No qualification of the data is necessary on DMCs RRF and %RSD/%D alone. Use professional judgment to evaluate DMCs and %RSD/%D data in conjunction with DMCs recoveries to determine the need for qualification of the data.

Qualify the initial calibration analytes listed in Table 2 using the following criteria in the CCVs:

Table 4. CCV Actions for Semivolatile Analysis

Criteria for Opening CCV	Criteria for Closing CCV	Action		
Criteria for Opening CCV	Criteria for Closing CCV	Detect	Non-detect	
CCV not performed at required frequency and sequence	CCV not performed at required frequency	Use professional judgment R	Use professional judgment R	
CCV not performed at specified concentration	CCV not performed at specified concentration			
RRF < Minimum RRF in Table 2 for target analyte	RRF < Minimum RRF in Table 2 for target analyte	Use professional judgment J or R	R	
RRF ≥ Minimum RRF in Table 2 for target analyte	RRF ≥ Minimum RRF in Table 2 for target analyte	No qualification	No qualification	
%D outside the Opening Maximum %D limits in Table 2 for target analyte	%D outside the Closing Maximum %D limits in Table 2 for target analyte	mits in Table 2 for target J		
%D within the inclusive Opening Maximum %D limits in Table 2 for target analyte	%D within the inclusive Closing Maximum %D limits in Table 2 for target analyte	No qualification	No qualification	

All criteria were met	_X	_
Criteria were not met		
and/or see below		_

BLANK ANALYSIS RESULTS (Sections 1 & 2)

The assessment of the blank analysis results is to determine the existence and magnitude of contamination problems. The criteria for evaluation of blanks apply only to blanks associated with the samples, including trip, equipment, and laboratory blanks. If problems with any blanks exist, all data associated with the case must be carefully evaluated to determine whether or not there is an inherent variability in the data for the case, or if the problem is an isolated occurrence not affecting other data.

List the contamination in the blanks below. High and low levels blanks must be treated separately.

Notes: The concentration of non-target compounds in all blanks must be less than or equal to 10 ug/L.

The concentration of target compounds in all blanks must be less than its CRQL listed in the method.

Samples taken from a drinking water tap do not have and associated field blank.

Laboratory blanks

Note:

DATE ANALYZED	LAB ID	LEVEL/ MATRIX	COMPOUND	CONCENTRATION UNITS
_No_target_ana	- llytes_detected_ -	_in_method_bla	anks	
Field/Equipmer	nt/Trip blank			
DATE ANALYZED	LAB ID	LEVEL/ MATRIX	COMPOUND	CONCENTRATION UNITS
_No_target_ana	lytes_detected_	_in_the_field/ed	uipment_blanks_analya	zed_with_this_data_package
				5, 420
	2000 de 1000 d			

12

All criteria were metX	_
Criteria were not met	
and/or see below	

BLANK ANALYSIS RESULTS (Section 3)

Blank Actions

Qualify samples based on the criteria summarized in Table 5:

Table 5. Blank and TCLP/SPLP LEB Actions for Semivolatile Analysis

Blank Type	Blank Result	Sample Result	Action
	Detect	Non-detect	No qualification
	< CRQL	< CRQL	Report at CRQL and qualify as non-detect (U)
		≥ CRQL	Use professional judgment
		< CRQL	Report at CRQL and qualify as non-detect (U)
Method,	≥CRQL	≥ CRQL but < Blank Result	Report at sample results and qualify as non-detect (U) or as unusable (R)
TCLP/SPLP LEB, Field		≥ CRQL and ≥ Blank Result	Use professional judgment
	Grossly high	Detect	Report at sample results and qualify as unusable (R)
	TIC > 5.0 ug/L (water) or 0.0050 mg/L (TCLP leachate) or TIC > 170 ug/Kg (soil)	Detect	Use professional judgment

List samples qualified

CONTAMINATION SOURCE/LEVEL	COMPOUND	CONC/UNITS	AL/UNITS	SQL	AFFECTED SAMPLES
					l
	,				

Note:

All criteria were metX_	
Criteria were not met	
and/or see below	

SURROGATE SPIKE RECOVERIES - DEUTERATED MONITORING COMPOUNDS (DMCs)

Laboratory performance of individual samples is established by evaluation of surrogate spike recoveries – deuterated monitoring compounds. All samples are spiked with surrogate compounds prior to sample analysis. The accuracy of the analysis is measured by the surrogate percent recovery. Since the effects of the sample matrix are frequently outside the control of the laboratory and may present relatively unique problems, the validation of data is frequently subjective and demands analytical experience and professional judgment.

Notes: Recoveries for DMCs in samples and blanks must be within the limits specified in Table 6

The recovery limits for any of the compounds listed in Table 6 may be expanded at any time during the period of performance if USEPA determines that the limits are too restrictive.

If a DMC is not added in the samples and blanks or the concentrations of DMCs in the samples and blank not the specified, use professional judgment in qualifying the data.

Table 7. DMC Actions for Semivolatile Analysis

Citaria	Action		
Criteria	Detect	Non-detect	
%R < 10% (excluding DMCs with 10% as a lower acceptance limit)	J-	R	
10% ≤ %R (excluding DMCs with 10% as a lower acceptance limit) < Lower Acceptance Limit	J-	υJ	
Lower Acceptance limit ≤ %R ≤ Upper Acceptance Limit	No qualification	No qualification	
%R > Upper Acceptance Limit	J+	No qualification	

SAMPLE ID	SURROGATE COMPOUND	ACTION
	uired_criteria_in_all_samples_analyzedNonde were_within_laboratory_recovery_limitsde	uterated_surrogates_added

Table 8. Semivolatile DMCs and the Associated Target Analytes

1,4-Dioxane-d ₈ (DMC-1)	Phenol-d ₅ (DMC-2)	Bis(2-Chloroethyl) ether-d _B (DMC-3)
1,4-Dioxane	Benzaldehyde	Bis(2-chloroethyl)ether
8	Phenol	2,2'-Oxybis(1-chloropropane)
		Bis(2-chloroethoxy)methane
2-Chlorophenol-d4(DMC-4)	4-Methylphenol-d ₈ (DMC-5)	4-Chloroaniline-d ₄ (DMC-6)
2-Chlorophenol	2-Methylphenol	4-Chloroaniline
	3-Methylphenol	Hexachlorocyclopentadiene
	4-Methylphenol	Dichlorobenzidine
	2,4-Dimethylphenol	
Nitrobenzene-d ₅ (DMC-7)	2-Nitrophenol-d ₄ (DMC-8)	2,4-Dichlorophenol-d3(DMC-9)
Acetophenone	Isophorone	2,4-Dichlorophenol
N-Nitroso-di-n-propylamine	2-Nitrophenol	Hexachlorobutadiene
l-lexachloroethane		Hexachlorocyclopentadiene
Nitrobenzene		4-Chloro-3-methylphenol
2,6-Dinitrotoluene		2,4,6-Trichlorophenol
2,4-Dinitrotoluene		2,4,5-Trichlorophenol
N-Nitrosodiphenylamine		1,2,4,5-Tetrachlorobenzene
		*Pentachlorophenol
		2,3,4,6-Tetrachlorophenol
Dimethylphthalate-d ₆ (DMC-10)	Acenaphthylene-da (DMC-11)	4-Nitrophenol-d4 (DMC-12)
Caprolactam	*Naphthalene	2-Nitroaniline
1,1'-Biphenyl	*2-Methylnaphthalene	3-Nitroaniline
Dimethylphthalate	2-Chloronaphthalene	2,4-Dinitrophenol
Diethylphthalate	*Acenaphthylene	4-Nitrophenol
Di-n-butylphthalate	*Acenaphthene	4-Nitroaniline
Butylbenzylphthalate		
Bis(2-ethylhexyl) phthalate		
Di-n-octylphthalate		

Fluorene-d ₁₀ (DMC-13)	4,6-Dinitro-2-methylphenol-d ₂ (DMC-14)	Anthracene-d ₁₀ (DMC-15)
Dibenzofuran *Fluorene	4,6-Dinitro-2-methylphenol	Hexachlorobenzene Atrazine
4-Chlorophenyl-phenylether		*Phenanthrene
4-Bromophenyl-phenylether		*Anthracene
Carbazole		
Pyrene-d ₁₀ (DMC-16)	Benzo(a)pyrene-d ₁₂ (DMC-17)	
*Fluoranthene	3,3'-Dichlorobenzidine	
*Pyrene	*Benzo(b)fluoranthene	
*Benzo(a)anthracene	*Benzo(k)fluoranthene	*
*Chrysene	*Benzo(a)pyrene	
	*Indeno(1,2,3-cd)pyrene	
	*Dibenzo(a,h)anthracene	
	*Benzo(g,h,i)perylene	

^{*}Included in optional Target Analyte List (TAL) of PAHs and PCP only.

Table 9. Semivolatile SIM DMCs and the Associated Target Analytes

Fluoranthene-d10 (DMC-1)	2-Methylnaphthalene-d10 (DMC-2)
Fluoranthene	Naphthalene
Pyrene	2-Methylnaphthalene
Benzo(a)anthracene	Acenaphthylene
Chrysene	Acenaphthene
Benzo(b)fluoranthene	Fluorene
Benzo(k)fluoranthene	Pentachlorophenol
Benzo(a)pyrene	Phenanthrene
Indeno(1,2,3-cd)pyrene	Anthracene
Dibenzo(a,h)anthracene	
Benzo(g,h,i)perylene	

All criteria were met _	
Criteria were not met	
and/or see below	_X

Method: SW846 8270D BY SIM

VII. A MATRIX SPIKE/MATRIX SPIKE DUPLICATE (MS/MSD)

This data is generated to determine long term precision and accuracy in the analytical method for various matrices. This data alone cannot be used to evaluate the precision and accuracy of individual samples. If any % R in the MS or MSD falls outside the designated range, the reviewer should determine if there are matrix effects, i.e. LCS data are within the QC limits but MS/MSD data are outside QC limit.

MS/MSD Recoveries and Precision Criteria

The laboratory should use one MS and a duplicate analysis of an unspiked field sample if target analytes are expected in the sample. If target analytes are not expected, MS/MSD should be analyzed.

NOTES:

Data for MS and MSDs will not be present unless requested by the Region. Notify the Contract Laboratory COR if a field or trip blank was used for the MS and MSD.

For a Matrix Spike that does not meet criteria, apply the action to only the field sample used to prepare the Matrix Spike sample. If it is clearly stated in the data validation materials that the samples were taken through incremental sampling or some other method guaranteeing the homogeneity of the sample group, then the entire sample group may be qualified.

List the %Rs, RPD of the compounds which do not meet the criteria.

Sample ID:	JC36372-4	Matrix/Level:	Soil
Sample ID:	JC36342-4_(SIM)	Matrix/Level:	Soil
Sample ID:	JC36371-1	Matrix/Level:	_Aqueous
Sample ID:	JC36371-1_(SIM)	Matrix/Level:	Aqueous
		91	E-11 12

The QC reported here applies to the following samples: JC36372-1, JC36372-2, JC36372-4

Compound Benzo(a)pyrene Dibenzo(a,h)	JC36373 ug/kg 135	2-4 Q	Spike ug/kg 37.9	MS ug/kg 220	MS % 224* a	Spike ug/kg 37.8	MSD ug/kg 564	MSD % 1134* a	RPD 88* b	Limits Rec/RPD 18-188/44
anthracene	23.2		37.9	59.9	97	37.8	118	251* a	65* b	28-169/41

⁽a) Outside of in house control limits due to possible sample nonhomogeneity.

Note: MS/MSD % recovery and RPD within laboratory control limits except for the cases described in this document. Results qualified as estimated (J or UJ) in sample JC36372-4 (SIM).

⁽b) Analytical precision exceeds in-house control limits.

^{* -} outside laboratory control limits

- * QC limits are laboratory in-house performance criteria, LL = lower limit, UL = upper limit.
- If QC limits are not available, use limits of 70 130 %.

Actions:

QUALITY	%R < LL	%R > UL
Positive results	J	J
Nondetects results	R	Accept

MS/MSD criteria apply only to the unspiked sample, its dilutions, and the associated MS/MSD samples:

If the % R for the affected compounds were < LL (or 70 %), qualify positive results (J) and nondetects (UJ).

If the % R for the affected compounds were > UL (or 130 %), only qualify positive results (J). If 25 % or more of all MS/MSD %R were < LL (or 70 %) or if two or more MS/MSD %Rs were < 10%, qualify all positive results (J) and reject nondetects (R).

A separate worksheet should be used for each MS/MSD pair.

All criteria were met __X__ Criteria were not met and/or see below ____

INTERNAL STANDARD PERFORMANCE

The assessment of the internal standard (IS) parameter is used to assist the data reviewer in determining the condition of the analytical instrumentation.

List the internal standard area of samples which do not meet the criteria.

DATE	SAMPLE ID	IS OUT	IS AREA	ACCEPTABLE RANGE	ACTION
Internal area	meets the requ	ired criteria for batch	samples corres	ponding to this data	package.

Action:

- 1. If an internal standard area count for a sample or blank is greater than 213.0% of the area for the associated standard (opening CCV or mid-point standard from initial calibration) (see Table 10 below):
 - a. Qualify detects for compounds quantitated using that internal standard as estimated low (J-).
 - b. Do not qualify non-detected associated compounds.
- 2. If an internal standard area count for a sample or blank is less than 20.0% of the area for the associated standard (opening CCV or mid-point standard from initial calibration):
 - a. Qualify detects for compounds quantitated using that internal standard as estimated high (J+).
 - b. Qualify non-detected associated compounds as unusable (R).
- 3. If an internal standard area count for a sample or blank is greater than or equal to 50.0%, and less than or equal to 213% of the area for the associated standard opening CCV or mid-point standard from initial calibration, no qualification of the data is necessary.
- 4. If an internal standard RT varies by more than 10.0 seconds: Examine the chromatographic profile for that sample to determine if any false positives or negatives exist. For shifts of a large magnitude, the reviewer may consider partial or total rejection of the data for that sample fraction. Detects should not need to be qualified as unusable (R) if the mass spectral criteria are met.
- 5. If an internal standard RT varies by less than or equal to 10.0 seconds, no qualification of the data is necessary.

Note: Inform the Contract Laboratory Program Project Officer (CLP PO) if the internal standard performance criteria are grossly exceeded. Note in the Data Review Narrative potential effects on the data resulting from unacceptable internal standard performance.

State in the Data Review Narrative if the required internal standard compounds are not added to a sample or blank or if the required internal standard compound is not analyzed at the specified concentration.

Actions:

Table 10. Internal Standard Actions for Semivolatile Analysis

Criteria -	Action		
CIRCIA	Detect	Non-detect	
Area response < 20% of the opening CCV or mid-point standard CS3 from ICAL	J÷	R	
20% ≤ Area response < 50% of the opening CCV or mid-point standard CS3 from ICAL	J+	נט	
50% ≤ Area response ≤ 200% of the opening CCV or mid-point standard CS3 from ICAL	No qualification	No qualification	
Area response > 200% of the opening CCV or mid-point standard CS3 from ICAL	J-	No qualification	
RT shift between sample/blank and opening CCV or mid-point standard CS3 from ICAL > 10.0 seconds	R	R	
RT shift between sample/blank and opening CCV or mid-point standard CS3 from ICAL < 10.0 seconds	No qualification	No qualification	

ang/or	r see below
TARGET COMPOUND IDENTIFICATION	
Criteria:	
Is the Relative Retention Times (RRTs) of reported compounds within ±0.06 I RRT [opening Continuing Calibration Verification (CCV) or mid-point s calibration].	
List compounds not meeting the criteria described above:	
Sample ID Compounds Action	
	
Mass spectra of the sample compound and a current laboratory-generated spectrum from the associated calibration standard (opening CCV or mid-p calibration)] must match according to the following criteria: a. All ions present in the standard mass spectrum at a relative in must be present in the sample spectrum. b. The relative intensities of these ions must agree within ±20% to sample spectra (e.g., for an ion with an abundance of 50% in the corresponding sample ion abundance must be between 30 c. lons present at greater than 10% in the sample mass spectrum standard spectrum, must be evaluated by a reviewer experimterpretation.	standard [i.e., the mass point standard from initial ntensity greater than 10% petween the standard and in the standard spectrum, 0-70%).
List compounds not meeting the criteria described above:	
Sample ID Compounds Action	ons
Identified compounds meet the required criteria	

Action:

- 1. The application of qualitative criteria for GC/MS analysis of target compounds requires professional judgment. It is up to the reviewer's discretion to obtain additional information from the laboratory. If it is determined that incorrect identifications were made, qualify all such data as unusable (R).
- 2. Use professional judgment to qualify the data if it is determined that cross-contamination has occurred.
- 3. Note in the Data Review Narrative any changes made to the reported compounds or concerns regarding target compound identifications. Note, for Contract Laboratory COR action, the necessity for numerous or significant changes.

TENTATIVELY IDENTIFIED COMPOUNDS (TICS)

NOTE: Tentatively identified compounds should only be evaluated when requested by a party from outside of the Hazardous Waste Support Section (HWSS).

List		Cs

Sample ID	Compound	Sample ID	Compound
=======================================			

Action:

- 1. Qualify all TIC results for which there is presumptive evidence of a match (e.g. greater than or equal to 85% match) as tentatively identified (NJ), with approximated concentrations. TICs labeled "unknown" are qualified as estimated (J).
- 2. General actions related to the review of TIC results are as follows:
 - a. If it is determined that a tentative identification of a non-target compound is unacceptable, change the tentative identification to "unknown" or another appropriate identification, and qualify the result as estimated (J).
 - b. If all contractually-required peaks were not library searched and quantitated, the Region's designated representative may request these data from the laboratory.
- 3. In deciding whether a library search result for a TIC represents a reasonable identification, use professional judgment. If there is more than one possible match, report the result as "either compound X or compound Y". If there is a lack of isomer specificity, change the TIC result to a nonspecific isomer result (e.g., 1,3,5-trimethyl benzene to trimethyl benzene isomer) or to a compound class (e.g., 2-methyl, 3-ethyl benzene to a substituted aromatic compound).
- 4. The reviewer may elect to report all similar compounds as a total (e.g., all alkanes may be summarized and reported as total hydrocarbons).

- 5. Target compounds from other fractions and suspected laboratory contaminants should be marked as "non-reportable".
- 6. Other Case factors may influence TIC judgments. If a sample TIC match is poor, but other samples have a TIC with a valid library match, similar RRT, and the same ions, infer identification information from the other sample TIC results.
- 7. Note in the Data Review Narrative any changes made to the reported data or any concerns regarding TIC identifications.
- 8. Note, for Contract Laboratory COR action, failure to properly evaluate and report TICs

All criteria were met _	X
Criteria were not met	
and/or see below	

SAMPLE QUANTITATION AND REPORTED CONTRACT REQUIRED QUANTITATION LIMITS (CRQLS)

Action:

- 1. When a sample is analyzed at more than one dilution, the lower CRQL are used unless a QC exceedance dictates the use of higher CRQLs from the diluted sample. Samples reported with an "E" qualifier should be reported from the diluted sample.
- 2. If any discrepancies are found, the Region's designated representative may contact the laboratory to obtain additional information that could resolve any differences. If a discrepancy remains unresolved, the reviewer must use professional judgment to decide which value is the most accurate. Under these circumstances, the reviewer may determine that qualification of data is warranted. Note in the Data Review Narrative a description of the reasons for data qualification and the qualification that is applied to the data.
- 3. For non-aqueous samples, if the solids is less than 10.0%, use professional judgment for both detects and non-detects. If the percent solid for a soil sample is greater than or equal to 10.0% and less than 30.0%, use professional judgment to qualify detects and non-detects. If the percent solid for a soil sample is greater than or equal to 30.0%, detects and non-detects should not be qualified (see Table 11).
- 4. Note, for Contract Laboratory COR action, numerous or significant failures to accurately quantify the target compounds or to properly evaluate and adjust CRQLs.
- 5. Results between MDL and CRQL should be qualified as estimated "J".
- 6. Results < MDL should be reported at the CRQL and qualified "U". MDLs themselves should not be reported.

Table 11. Percent Solids Actions for Semivolatile Analysis for Non-Aqueous Samples

Criteria	Ac	Action			
Cincia	Detects	Non-detects			
%Solids < 10.0%	Use professional judgment	Use professional judgment			
10.0% ≤ %Solids ≤ 30.0%	Use professional judgment	Use professional judgment			
%Solids > 30.0%	No qualification	No qualification			

SAMPLE QUANTITATION

The sample quantitation evaluation is to verify laboratory quantitation results. In the space below, please show a minimum of one sample calculation:

Sample ID:_	_ JC363	72-1	Analyte:Chrysene	RF:_1.185_
[]	=	(491575)(40)	/(1034286)(1.185)	
	=	16.04 ppm	Ok	

QUANTITATION LIMITS

A. Dilution performed

SAMPLE ID	DILUTION FACTOR	REASON FOR DILUTION
	-	
2		

				Criter	teria were metN/A ia were not met r see below	
FIELD DUPLICATE	PRECIS	SION				
Sample IDs				Mat	rix:	_
analyses measure to laboratory duplicates will have a greater field duplicate sample. The project QAPP structure of the structure of the project content of the	both field s which variance les. hould be if large I	d and lab precision only laboratory per than water matrice reviewed for project RPD (> 50 %) is a second or project than the control of the c	en; therefore, the result erformance. It is also exces due to difficulties and exceptions.	Its may hexpected associated.	overall precision. The nave more variability the that soil duplicate resed with collecting ident of the samples and rebled.	nan ults ical
COMPOUND	SQL ug/L	SAMPLE CONC. (ug/l)	DUPLICATE CONC. (ug/l)	RPD	ACTION	
used to assess prec	ision. RI	PD within the requi	of this data package ired guidance docume ses described in this do	nt criteria	•	

All criteria were met _	_X_	
Criteria were not met		
and/or see below	22	

OTHER ISSUES

ziot oampioo quami	ed based on the degradation of system	r performance during simple analysis.
Sample ID	Comments	Actions
		· ·
Action:		
during sample and degradation of syst	alyses. Inform the Contract Laborato em performance which significantly aff	mined that system performance has degraded by Program COR any action as a result of ected the data.
B. Overall Ass	essment of Data	
List samples qualifi	ed based on other issues:	
Sample ID	Comments	Actions
		e_dataResults_are_valid_and_can_be_used wn_below

Action:

- 1. Use professional judgment to determine if there is any need to qualify data which were not qualified based on the Quality Control (QC) criteria previously discussed.
- Write a brief narrative to give the user an indication of the analytical limitations of the data. Inform the Contract Laboratory COR the action, any inconsistency of the data with the Sample Delivery Group (SDG) Narrative. If sufficient information on the intended use and required quality of the data is available, the reviewer should include their assessment of the usability of

- the data within the given context. This may be used as part of a formal Data Quality Assessment (DQA).
- 3. Sometimes, due to dilutions, re-analysis or SIM/Scan runs are being performed, there will be multiple results for a single analyte from a single sample. The following criteria and professional judgment are used to determine which result should be reported:
 - The analysis with the lower CRQL
 - The analysis with the better QC results
 - The analysis with the higher results

EXECUTIVE NARRATIVE

SDG No:

JC36372

Laboratory:

Accutest, New Jersey

Analysis:

SW846-8015C

Number of Samples:

Ω

Location:

BMSMC, Humacao, PR

SUMMARY:

Eight (8) samples were analyzed for the low molecular weight alcohols (LMWAs) list following method SW846-8015C. The sample results were assessed according to USEPA data validation guidance documents in the following order of precedence: "Test Methods for Evaluating Solid Waste, Physical/Chemical Methods SW-846 (Final Update III, December 1996)," specifically for Methods 8000/8015C are utilized. The QC criteria and data validation actions listed on the data review worksheets are from the primary guidance document, unless otherwise noted.

Results are valid and can be used for decision making purposes.

Critical issues:

None

Major:

None

Minor:

None

Critical findings:

None

Major findings:

None

Minor findings:

1. MS/MSD % recovery and RPD within laboratory control limits except for the cases described in the Data Review Worksheet. RPD for ethanol, isopropyl alcohol, and methanol outside the laboratory control limits. No qualification made based on RPD results, professional judgment.

No MS/MSD samples analyzed for the aqueous matrix. BS/BSD used to assess accuracy.

COMMENTS:

Results are valid and can be used for decision making purposes.

Reviewers Name:

Rafael Infante

Chemist License 1888

Rafael Infact

Signature:

Date:

February 18, 2017

SAMPLE ORGANIC DATA SAMPLE SUMMARY

Sample ID: JC36372-1

Sample location: BMSMC, Humacao, PR

Sampling date: 1/26/2017

Matrix: Soil

METHOD: 8015C

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Ethanol	120	ug/kg	1.0	-	U	Yes
Isobutyl Alcohol	120	ug/kg	1.0	-	U	Yes
Isopropyl Alcohol	120	ug/kg	1.0	•	U	Yes
n-Propyl Alcohol	120	ug/kg	1.0	-	U	Yes
n-Butyl Alcohol	120	ug/kg	1.0	-	U	Yes
sec-Butyl Alcohol	120	ug/kg	1.0	-	U	Yes
Methanol	230	ug/kg	1.0	•	U	Yes

Sample ID: JC36372-2

Sample location: BMSMC, Humacao, PR

Sampling date: 1/26/2017

Matrix: Soil

METHOD: 8015C

WILTHOU.	,0250					
Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Ethanol	120	ug/kg	1.0	-	U	Yes
Isobutyl Alcohol	120	ug/kg	1.0	•	U	Yes
Isopropyl Alcohol	120	ug/kg	1.0	-	U	Yes
n-Propyl Alcohol	120	ug/kg	1.0	-	U	Yes
n-Butyl Alcohol	120	ug/kg	1.0	-	U	Yes
sec-Butyl Alcohol	120	ug/kg	1.0	-	U	Yes
Methanol	240	ug/kg	1.0	-	U	Yes

Sample ID: JC36372-3

Sample location: BMSMC, Humacao, PR

Sampling date: 1/26/2017

Matrix: Soil

METHOD: 8015C

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Ethanol	120	ug/kg	1.0	•	U	Yes
Isobutyl Alcohol	120	ug/kg	1.0	-	U	Yes
Isopropyl Alcohol	120	ug/kg	1.0	-	U	Yes
n-Propyl Alcohol	120	ug/kg	1.0	-	υ	Yes
n-Butyl Alcohol	120	ug/kg	1.0	•	U	Yes
sec-Butyl Alcohol	120	ug/kg	1.0	-	U	Yes
Methanol	250	ug/kg	1.0	_	U	Yes

Sample ID: JC36372-4

Sample location: BMSMC, Humacao, PR

Sampling date: 1/26/2017

Matrix: Soil

METHOD: 8015C

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Ethanol	120	ug/kg	1.0	o 7 .0	U	Yes
Isobutyl Alcohol	120	ug/kg	1.0	-	U	Yes
Isopropyl Alcohol	120	ug/kg	1.0	-	U	Yes
n-Propyl Alcohol	120	ug/kg	1.0		U	Yes
n-Butyl Alcohol	120	ug/kg	1.0	-	U	Yes
sec-Butyl Alcohol	120	ug/kg	1.0	-	U	Yes
Methanol	230	ug/kg	1.0	-	U	Yes

Sample 1D: JC36372-5

Sample location: BMSMC, Humacao, PR

Sampling date: 1/26/2017

Matrix: AQ - Field Blank Soil

METHOD: 8015C

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Ethanol	100	ug/l	1.0	-	U	Yes
Isobutyl Alcohol	100	ug/l	1.0	-	U	Yes
Isopropyl Alcohol	100	ug/l	1.0	-	U	Yes
n-Propyl Alcohol	100	ug/l	1.0	-	U	Yes
n-Butyl Alcohol	100	ug/l	1.0	-	U	Yes
sec-Butyl Alcohol	100	ug/l	1.0	-	U	Yes
Methanol	200	ug/l	1.0	-	U	Yes

Sample ID: JC36372-6

Sample location: BMSMC, Humacao, PR

Sampling date: 1/26/2017

Matrix: AQ - Equipment Blank

METHOD: 8015C

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Ethanol	100	ug/l	1.0	-	U	Yes
Isobutyl Alcohol	100	ug/l	1.0	-	U	Yes
Isopropyl Alcohol	100	ug/l	1.0	-	U	Yes
n-Propyl Alcohol	100	ug/l	1.0	-	U	Yes
n-Butyl Alcohol	100	ug/l	1.0	-	U	Yes
sec-Butyl Alcohol	100	ug/l	1.0	-	U	Yes
Methanol	200	ug/l	1.0	-	U	Yes

Sample ID: JC36372-4MS

Sample location: BMSMC, Humacao, PR

Sampling date: 1/26/2017

Matrix: Soil

METHOD: 8015C

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Ethanol	6290	ug/kg	1.0	-	-	Yes
Isobutyl Alcohol	6650	ug/kg	1.0	-	-	Yes
Isopropyl Alcohol	6580	ug/kg	1.0	-	+3	Yes
n-Propyl Alcohol	6490	ug/kg	1.0	-	-	Yes
n-Butyl Alcohol	6730	ug/kg	1.0	-	7.5	Yes
sec-Butyl Alcohol	7450	ug/kg	1.0	-	20	Yes
Methanol	6430	ug/kg	1.0	-	7.0	Yes

Sample ID: JC36372-4MSD

Sample location: BMSMC, Humacao, PR

Sampling date: 1/26/2017

Matrix: Soil

METHOD: 8015C

METHOD:	POTOC					
Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Ethanol	4730	ug/kg	1.0	-	-	Yes
Isobutyl Alcohol	5550	ug/kg	1.0	-	-	Yes
Isopropyl Alcohol	4730	ug/kg	1.0	-	-	Yes
n-Propyl Alcohol	7660	ug/kg	1.0	-		Yes
n-Butyl Alcohol	6190	ug/kg	1.0		-	Yes
sec-Butyl Alcohol	6690	ug/kg	1.0	-	-	Yes
Methanol	3890	ug/kg	1.0	-	*	Yes

	Project Number:JC36372
	Date:01/26/2017
	Shipping Date:01/26/2017
	EPA Region: 2
	<u></u>
REVIEW OF VOLATILE O The following guidelines for evaluating volatile organics were document will assist the reviewer in using professional judg serving the needs of the data users. The sample results guidance documents in the following order of precede Physical/Chemical Methods SW-846 (Final Update III, Decer utilized. The QC criteria and data validation actions listed guidance document, unless otherwise noted. The hardcopied (laboratory name) _Accutest and the quality control and performance data summarized. The	e created to delineate required validation actions. This gment to make more informed decision and in bette were assessed according to USEPA data validation ince: "Test Methods for Evaluating Solid Wasternber 1996)," specifically for Methods 8000/8015C are on the data review worksheets are from the primary data package received has been reviewed
Lab. Project/SDG No.:JC36372	Sample matrix:Soil
Trip blank No.:	
Overall Comments:_Low_molecular_weight_alcohols_b	y_SW-846_8015C
Definition of Qualifiers: J- Estimated results U- Compound not detected R- Rejected data UJ- Estimated nondetect Reviewer: Date:February_18,_2017	

DATA COMPLETENESS

MISSING INFORMATION	DATE LAB. CONTACTED	DATE RECEIVED
	4 M - 4 A A A	
	141	
	<u> </u>	
·		
		1
10,		
-		_

All criteria were metX	
Criteria were not met	
and/or see below	

HOLDING TIMES

The objective of this parameter is to ascertain the validity of the results based on the holding time of the sample from time of collection to the time of analysis.

Complete table for all samples and note the analysis and/or preservation not within criteria

SAMPLE ID	DATE SAMPLED	DATE ANALYZED	pΗ	ACTION
All samples analy	zed within the recomr	nended method holding	. All sam	ples properly preserved.
			1	
			1	4

<u>Criteria</u>

Aqueous samples – 14 days from sample collection for preserved samples (pH \leq 2, 4°C), no air bubbles. Aqueous samples – 7 days from sample collection for unpreserved samples, 4°C, no air bubbles.

Soil samples- 7 days from sample collection.

Cooler temperature (Criteria: 4 ± 2 °C): 3.6°C

Actions

If the VOCs vial(s) have air bubbles, estimate positive results (J) and reject nondetects (R).

If the % solids of soil samples is 10-50%, estimates positive results (J) and nondetects (UJ)

If the % solid of soil samples is < 10%, estimate positive results (J) and reject nondetects (R).

If holding times are exceeded but < 14 days beyond criteria, estimate positive results (J) and nondetects (UJ).

If holding times are exceeded but < 28 days beyond criteria, estimate positive results (J) and reject nondetects (R).

If holding times are grossly exceeded (> 28 days beyond criteria), reject all results (R).

If samples were not iced or if the ice were melted (> 10°C), estimate positive results (J) and nondetects (UJ).

	All criteria were metN/A	
Criteria	were not met see helow	

GC/MS TUNING

The assessment of the tuning results is to determine if the sample instrumentation is within the standard tuning QC limits
N/A_ The BFB performance results were reviewed and found to be within the specified criteria.
N/A_ BFB tuning was performed for every 12 hours of sample analysis.
If no, use professional judgment to determine whether the associated data should be accepted, qualified or rejected.
List the samples affected:
If mass calibration is in error, all associated data are rejected.

All criteria were met _	х
Criteria were not met	
and/or see below	

CALIBRATION VERIFICATION

Compliance requirements for satisfactory instrument calibration are established to ensure that the instrument is capable of producing and maintaining acceptable quantitative data.

	Dat	e of initial calibration:	10/10/16		
	Dat	es of continuing calibra	tion:01/30/17		
	Dat	es of final calibration ve	erification:10/10/1	0;_01/30/17	
	Inst	trument ID number:	GCGH_		
	Mai	trix/Level:	Aqueous/low		
DATE	LAB FILE ID#	CRITERIA OUT RFs, %RSD, %D, r	COMPOUND	SAMPLES AFFECTED	

Note:	Initial, continuing, and final calibration verifications meets method specific criteria in at least one
	of the two columns. Final calibration verification included in data packages.

Criteria

All RFs must be > 0.05 regardless of method requirements for SPCC.

All %RSD must be \leq 15 % regardless of method requirements for CCC.

All %Ds must be ≤ 20% regardless of method requirements for CCC.

It should be noted that Region 2 SOP HW-24 does not specify criterion for the curve correlation coefficient (r). A limit for r of \geq 0.995 has therefore been utilized as professional judgment.

Actions

If any compound has an initial RF or a continuing RF of < 0.05, estimate positive results (J) and reject nondetects (R), regardless of method requirements.

If any compound has a %RSD > 15%, estimate positive results (J) and use professional judgment to qualify nondetects.

If any compound has a %RSD > 90%, estimate positive results (J) and reject nondetects (R).

If any compound has a % D > 20%, estimate positive results (J) and reject nondetects (R).

If any compound has a % D > 20%, estimate positive results (J) and nondetects (UJ).

If any compound has a % D > 90%, estimate positive results (J) and reject nondetects (R).

If any compound has r < 0.995, estimate positive results and nondetects.

A separate worksheet should be filled for each initial curve

All criteria were metX	
Criteria were not met	
and/or see below	

V A. BLANK ANALYSIS RESULTS (Sections 1 & 2)

The assessment of the blank analysis results is to determine the existence and magnitude of contamination problems. The criteria for evaluation of blanks apply only to blanks associated with the samples, including trip, equipment, and laboratory blanks. If problems with any blanks exist, all data associated with the case must be carefully evaluated to determine whether or not there is an inherent variability in the data for the case, or if the problem is an isolated occurrence not affecting other data.

List the contamination in the blanks below. High and low levels blanks must be treated separately.

Laboratory blanks

DATE ANALYZED	LAB ID	LEVEL/ MATRIX	COMPOUND	CONCENTRATION UNITS
			fic_criteria	
Field/Equipme	<u>nt</u> /Trip blank			
DATE ANALYZED	LAB ID	LEVEL/ Matrix	COMPOUND	CONCENTRATION UNITS
				zed_with_this_data_package
- 350			V	

All criteria were met _	_X
Criteria were not met	
and/or see below	

VB. BLANK ANALYSIS RESULTS (Section 3)

Blank Actions

Action Levels (ALs) should be based upon the highest concentration of contaminant determined in any blank. Do not qualify any blank with another blank. The ALs for samples which have been diluted should be corrected for the sample dilution factor and/or % moisture, where applicable. No positive sample results should be reported unless the concentration of the compound in the samples exceeds the ALs:

ALs = 10x the amount of common contaminants (methylene chloride, acetone, 2-butanone, and toluene) ALs = 5x for any other compounds

Specific actions are as follows:

If the concentration is < sample quantitation limit (SQL) and \le AL, report the compound as not detected (U) at the SQL.

If the concentration is \geq SQL but \leq AL, report the compound as not detected (U) at the reported concentration.

If the concentration is ≥ SQL and > AL, report the concentration unqualified.

Notes:

High and low level blanks must be treated separately Compounds qualified "U" for blank contamination are still considered "hits" when qualifying for calibration criteria.

CONTAMINATION SOURCE/LEVEL	COMPOUND	CONC/UNITS	AL/UNITS	SQL	AFFECTED SAMPLES
				- T	
	-				

All criteria were met	_X	
Criteria were not met		
and/or see below		

SURROGATE SPIKE RECOVERIES

Laboratory performance of individual samples is established by evaluation of surrogate spike recoveries. All samples are spiked with surrogate compounds prior to sample analysis. The accuracy of the analysis is measured by the surrogate percent recovery. Since the effects of the sample matrix are frequently outside the control of the laboratory and may present relatively unique problems, the validation of data is frequently subjective and demands analytical experience and professional judgment.

List the percent recoveries (%Rs) which do not meet the criteria for surrogate recovery.

Matrix: solid/aqueous

SAMPLE ID		SL	JRROGATE	COMPOUND		ACTION
	Hexan	ol DB I	FM	TOL-d8	BFB	
	S1 a	S1 b				
JC36372-5	111	94				
JC36372-6	115	102				
GGH5639-BS	102	102				
GGH5639-BSD	102	105				
GGH5639-MB1	99	111				
JC36372-1	95	93				
JC36372-2	85	87				
JC36372-3	101	104				
JC36372-4	76	71				
GGH5640-BS	100	94				
GGH5640-MB1	96	92				
JC36372-4MS	112	109				
JC36372-4MSD	93	94				
(a) Recovery from GC s	ignal #2		(b) Recove	ery from GC sign	al #1	
Note: All surre	ogate re	coveries	within labo	ratory control lim	its.	

QC Limits* (Aqueous)				
LL_to_UL	_56_to_145_	to	to	to
QC Limits* (Solid-Low)				
LL_to_UL	_52_to_141	to	to	to
QC Limits* (Solid-Med)				
LL_to_UL	to	to	to	to
1,2-DCA = $1,2$ -Dichloro	methane-d4		TOL-d8 = T	oluene-d8
DBFM = Dibromofluoro	methane		BFB = Bror	mofluorobenzene

- * QC limits are laboratory in-house performance criteria, LL = lower limit, UL = upper limit.
- * If QC limits are not available, use limits of 80 120 % for aqueous and 70 130 % for solid samples.

Actions:

QUALITY	%R < 10%	%R = 10% - LL	%R > UL
Positive results	J	J	J
Nondetects results	R	UJ	Accept

Surrogate action should be applied:

If one or more surrogate in the VOC fraction is out of specification, but has a recovery of > 10%. If any one surrogate in a fraction shows < 10 % recovery.

All criteria were metX
Criteria were not met
and/or see below

VII. A MATRIX SPIKE/MATRIX SPIKE DUPLICATE (MS/MSD)

This data is generated to determine long term precision and accuracy in the analytical method for various matrices. This data alone cannot be used to evaluate the precision and accuracy of individual samples. If any % R in the MS or MSD falls outside the designated range, the reviewer should determine if there are matrix effects, i.e. LCS data are within the QC limits but MS/MSD data are outside QC limit.

1. MS/MSD Recoveries and Precision Criteria

The laboratory should use one MS and a duplicate analysis of an unspiked field sample if target analytes are expected in the sample. If target analytes are not expected, MS/MSD should be analyzed.

List the %Rs, RPD of the compounds which do not meet the criteria.

Sample ID:JC36372-4MS/-4MSD				Matrix/Level:Soil/low		
MS OR MSD	COMPOUND	% R	RPD	QC LIMITS	ACTION	
	ecoveries_and_RPD_ent				cept_for_the_cases_d	escribed_

Note: RPD for ethanol, isopropyl alcohol, and methanol outside the laboratory control limits. No qualification made based on RPD results, professional judgment.

No MS/MSD samples analyzed for the aqueous matrix. BS/BSD used to assess accuracy.

- * QC limits are laboratory in-house performance criteria, LL = lower limit, UL = upper limit.
- * If QC limits are not available, use limits of 70 130 %.

Actions:

QUALITY	%R < LL	%R > UL
Positive results	J	J
Nondetects results	R	Accept

All criteria were met __X__ Criteria were not met and/or see below ____

MS/MSD criteria apply only to the unspiked sample, its dilutions, and the associated MS/MSD samples:

If the % R for the affected compounds were < LL (or 70 %), qualify positive results (J) and nondetects (UJ).

If the % R for the affected compounds were > UL (or 130 %), only qualify positive results (J). If 25 % or more of all MS/MSD %R were < LL (or 70 %) or if two or more MS/MSD %Rs were < 10%, qualify all positive results (J) and reject nondetects (R).

VII. B MATRIX SPIKE/MATRIX SPIKE DUPLICATE

MS/MSD - Unspiked Compounds

It should be noted that Region 2 SOP HW-24 does not specify a MS/MSD criteria for the unspiked compounds in the sample. A %RSD of < 50% has therefore been utilized as professional judgment.

If all target analytes were spiked in the MS/MSD, this review element is not applicable.

List the %RSD of the compounds which do not meet the criteria.

Sample ID:	-		Matrix/Le	evel/Unit:		
COMPOUND	SAMPLE CONC.	MS CONC.	MSD CONC.	% RSD	ACTION	
	100		-3000			
Table 1						

Actions:

A separate worksheet should be used for each MS/MSD pair.

^{*} If the % RSD > 50, qualify the positive result in the unspiked samples as estimated (J).

^{*} If the % RSD is not calculated (NC) due to nondetected value, use professional judgment to qualify the data.

All criteria were met _	_X	
Criteria were not met		
and/or see below		

VIII. LABORATORY CONTROL SAMPLE (LCS) ANALYSIS

This data is generated to determine accuracy of the analytical method for various matrices.

1. LCS Recoveries Criteria

Where LCS spiked with the same analyte at the same concentrations as the MS/MSD? Yes or No. If no make note in data review memo.

List the %R of compounds which do not meet the criteria

_	

Note:

- * QC limits are laboratory in-house performance criteria, LL = lower limit, UL = upper limit.
- * If QC limits are not available, use limits of 70 130 %.

Actions:

QUALITY	%R < LL	%R > UL
Positive results	J	J
Nondetects results	R	Accept

All analytes in the associated sample results are qualified for the following criteria.

If 25 % of the LCS recoveries were < LL (or 70 %), qualify all positive results (j) and reject nondetects (R).

If two or more LCS were below 10 %, qualify all positive results as (J) and reject nondetects (R).

2. Frequency Criteria:

Where LCS analyzed at the required frequency and for each matrix? <u>Yes</u> or No. If no, the data may be affected. Use professional judgment to determine the severity of the effect and qualify data accordingly. Discuss any actions below and list the samples affected.

All criteria were met	_X_
Criteria were not met	
and/or see below	

IX. FIELD/LABORATORY DUPLICATE PRECISION

Sample IDs:JC36372-2/JC36372-3	Matrix:Soil
--------------------------------	-------------

Field/laboratory duplicates samples may be taken and analyzed as an indication of overall precision. These analyses measure both field and lab precision; therefore, the results may have more variability than laboratory duplicates which only laboratory performance. It is also expected that soil duplicate results will have a greater variance than water matrices due to difficulties associated with collecting identical field duplicate samples.

The project QAPP should be reviewed for project-specific information.

Suggested criteria: RPD \pm 30% for aqueous samples, RPD \pm 50 % for solid samples. If both samples and duplicate are <5 SQL, the RPD criteria is doubled.

COMPOUND	SQL	SAMPLE CONC.	DUPLICATE CONC.	RPD	ACTION
Field duplicate	e analys	rad with this data no	okago PPD within Joha	rotoni a	enerally acceptable and
l leid duplicate			performance criteria cor		
					·

Actions:

Qualify as estimated positive results (J) and nondetects (UJ) for the compound that exceeded the above criteria. For organics, only the sample and duplicate will be qualified.

If an RPD cannot be calculated because one or both of the sample results is not detected, the following actions apply:

If one sample result is not detected and the other is greater than 5x the SQL qualify (J/UJ).

If one sample value is not detected and the other is greater than 5x the SQL and the SQLs for the sample and duplicate are significantly different, use professional judgment to determine if qualification is appropriate.

If one sample value is not detected and the other is less than 5x, use professional judgment to determine if qualification is appropriate.

If both sample and duplicate results are not detected, no action is needed.

All criteria were metN/	Α
Criteria were not met	
and/or see below	

X. INTERNAL STANDARD PERFORMANCE

The assessment of the internal standard (IS) parameter is used to assist the data reviewer in determining the condition of the analytical instrumentation.

List the internal standard area of samples which do not meet the criteria.

- * Area of +100% or -50% of the IS area in the associated calibration standard.
- * Retention time (RT) within 30 seconds of the IS area in the associated calibration standard.

DATE	SAMPLE ID	IS OUT	IS AREA	ACCEPTABLE ACTION RANGE	
					121
		- 5	900		
	TO SERVICE STATE OF THE PERSON NAMED IN COLUMN TO SERVICE STATE OF THE PERSON NAMED				
- 5					

Actions:

1. IS actions should be applied to the compound quantitated with the out-of-control ISs

QUALITY	IS AREA < -25%	IS AREA = -25 % TO - 50%	IS AREA > + 100%
Positive results	J	J	J
Nondetected results	R	UJ	ACCEPT

2. If a IS retention time varies more than 30 seconds, the chromatographic profile for that sample must be examined to determine if any false positive or negative exists. For shifts of a large magnitude, the reviewer may consider partial or total rejection of the data for the sample fraction.

All criteria were met _	_X	_
Criteria were not met		
and/or see below		

XII. SAMPLE QUANTITATION

The sample quantitation evaluation is to verify laboratory quantitation results. In the space below, please show a minimum of one sample calculation:

JC36372-4MS

Sec-Butyl alcohol

RF = 26.42

All criteria were met	Χ_	_
Criteria were not met		
and/or see below		

XII. QUANTITATION LIMITS

A. Dilution performed

SAMPLE ID	DILUTION FACTOR	REASON FOR DILUTION
		Y t
		1/2
		-
	7/2	
	7,5	

B. Pe	cent Solids
-------	-------------

List sam	ples whic	h have <u>≤</u> \$	50 % solid	ds			
			 -			1555	9
Table 1							

Actions:

If the % solids of a soil sample is 10-50%, estimate positive results (J) and nondetects (UJ)

If the % solids of a soil sample is < 10%, estimate positive results (J) and reject nondetects (R)

EXECUTIVE NARRATIVE

SDG No:

JC36372

Laboratory:

Accutest, New Jersey

Analysis:

SW846-8081B

Number of Samples:

8

Location:

BMSMC, Humacao, PR

SUMMARY:

Eight (8) samples were analyzed for the TCL pesticides list following method SW846-8081B. The sample results were assessed according to USEPA data validation guidance documents in the following order of precedence *Hazardous Waste Support Section SOP No. HW-36A, Revision O, June, 2015. SOM02.2. Pesticide Data Validation.* The QC criteria and data validation actions listed on the data review worksheets are from the primary guidance document upless otherwise noted.

guidance document, unless otherwise noted.

Results are valid and can be used for decision making purposes.

Critical issues:

None

Major:

None

Minor:

None

Critical findings:

None

Major findings:

None

Minor findings:

- 1. Initial and initial calibration verification within the guidance document performance criteria. Continuing calibration % differences meet the performance criteria in at least one of the two columns. Final calibration verification not included in data package. No action taken.
- 2. MS/MSD % recoveries and RPD within laboratory control limits except for the cases described in the Data Review Worksheet. Results for Aldrin qualified as estimated (J) in sample JC36372-4. No qualifications made based on RPD results, professional judgment.

COMMENTS:

Results are valid and can be used for decision making purposes.

Reviewers Name:

Rafael Infante

Chemist License 1888

Signature:

Date:

February 18, 2017

SAMPLE ORGANIC DATA SAMPLE SUMMARY

Sample ID: JC36372-1
Sample location: BMSMC, Humacao, PR
Sampling date: 26-Jan-17
Matrix: Soil

	Reportable	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
	Lab Flag Validation	n	n	n	n	n	n	n	n	n	n	n	n	D	n	-	n	⊃	⊃	⊃	⊃	n
	Lab Flag	1	ı	ı	•	1	ı	ı	1		ı	1	ı	1	•	•	•	•		•	1	i
	Jnits Dilution Factor	1	1	1	1	П	1	1	1	1	П	П	П	ਜ	П	₽	П	₽	+	1	1	Т
	Units [ug/kg	ug/kg	ug/kg	ug/kg	ug/kg	ug/kg	ug/kg	ug/kg	ug/kg	ug/kg	ug/kg	ug/kg	ug/kg	ug/kg	ug/kg	ug/kg	ug/kg	ug/kg	ug/kg	ug/kg	ug/kg
VEHICLE COOLE	Result	0.76	0.76	0.76	0.76	0.76	0.76	0.76	0.76	0.76	0.76	0.76	0.76	0.76	0.76	0.76	0.76	0.76	0.76	1.5	0.76	19
-	Analyte Name	Aldrin	alpha-BHC	beta-BHC	delta-BHC	gamma-BHC (Lindane)	alpha-Chlordane	gamma-Chlordane	Dieldrin	4,4'-DDD	4,4'-DDE	4,4'-DDT	Endrin	Endosulfan sulfate	Endrin aldehyde	Endosulfan-1	Endosulfan-II	Heptachlor	Heptachlor epoxide	Methoxychlor	Endrin ketone	Toxaphene

Sample ID: JC36372-2 Sample location: BMSMC, Humacao, PR

Sampling date: 26-Jan-17 Matrix: Soil

	Reportable	Yes	Yes	Yes	Yes		Yes	Yes	Yes Yes	Yes Yes Yes	Yes Yes Yes	Yes Yes Yes	Yes Yes Yes Yes	Yes Yes Yes Yes	Yes Yes Yes Yes Yes	Yes Yes Yes Yes	Yes Yes Yes Yes Yes	Yes Yes Yes Yes Yes	Yes	Yes Yes Yes Yes Yes Yes Yes Yes	Yes Yes Yes Yes Yes Yes Yes Yes Yes	Yes
	Lab Flag Validation	⊃	ח	ם	ם	כ	כ	Þ	⊃	כ	⊃)	コ	כ	כ	o	ם	ם	כ	n	n	
	Lab Flag	,		•	•	â		1	ι	,	1	1	•	•	•	•	ī	•	,	70	×	
	Jnits Dilution Factor	П	П	Н	П	П	П	П	П	П	П	П	П	П	П	П		П	П	П	н	•
	Units	ug/kg	ug/kg	ug/kg	ug/kg	ug/kg	ug/kg	ug/kg	ug/kg	ug/kg	ug/kg	ug/kg	ug/kg	ug/kg	A.m.							
	#					_	_	_	_						_			_	-+		_	
METHOD, SUSTE	Result	0.84	0.84	0.84	0.84	0.84	0.84	0.8	0.8	0.84	0.84	0.84	0.84	0.84	0.84	0.84	0.84	78.0	0.8	1.7	0.8	1,0

Sample ID: JC36372-3
Sample location: BMSMC, Humacao, PR
Sampling date: 26-Jan-17
Matrix: Soil

1	Units Dilution Factor
1	ug/kg
1 - 0 Yes	ug/kg
1 - 0 Yes	ug/kg
1 - 0 Yes 1 1 - 0 Yes	ug/kg
1	ug/kg
1	ug/kg
1 - 0 Yes	ug/kg
1 - U Yes 1 - U	ug/kg
1	ug/kg
1 - U Yes 1	ug/kg
1 Yes 1 Yes 1 Yes 1 Yes 1	ug/kg
1 Yes U Yes 1	ug/kg
1 Ves	ug/kg
	ug/kg

Sample ID: JC36372-4

Sample location: BMSMC, Humacao, PR

Sampling date: 26-Jan-17 Matrix: Soil

	1 1 0 0 0 0						
Analyte Name	Result	Units Dilu	Jnits Dilution Factor	Lab Flag	Lab Flag Validation	Reportable	
Aldrin	0.76	ug/kg	-	,	m	Yes	
alpha-BHC	0.76	ug/kg	-	1	n	Yes	
beta-BHC	0.76	ug/kg	1	1	⊃	Yes	
delta-BHC	0.76	ug/kg	-	1	n	Yes	
gamma-BHC (Lindane)	0.76	ug/kg	1	6	⊃	Yes	
alpha-Chlordane	0.76	ug/kg	1	1	⊃	Yes	
gamma-Chlordane	0.76	ug/kg	1	1	⊃	Yes	
Dieldrin	0.76	ug/kg	_	•	⊃	Yes	
4,4'-DDD	0.76	ug/kg		,	n	Yes	
4,4'-DDE	0.76	ug/kg	Т	1	⊃	Yes	
4,4'-DDT	0.76	ug/kg	Т	•	⊃	Yes	
Endrin	0.76	ug/kg	₽	,	⊃	Yes	
Endosulfan sulfate	0.76	ug/kg	₽	1	⊃	Yes	
Endrin aldehyde	0.76	ug/kg	₽	1	⊃	Yes	
Endosulfan-I	0.76	ug/kg	н	٠	⊃	Yes	
Endosulfan-11	0.76	ug/kg	↔	,	5	Yes	
Heptachlor	0.76	ug/kg	Н	r	⊃	Yes	
Heptachlor epoxide	0.76	ug/kg	П	٠	⊃	Yes	
Methoxychlor	1.5	ug/kg	Н	•	ח	Yes	
Endrin ketone	0.76	ug/kg	Н	•	D	Yes	
Toxaphene	19	ug/kg	—	ř.	n	Yes	

Sample ID: JC36372-5

Sample location: BMSMC, Humacao, PR Sampling date: 26-Jan-17 Matrix: AQ - Field Blank Soil

Analyte Name	Result	Units Dilution Factor	on Factor	Lab Flag	Lab Flag Validation	Reportable
Aldrin	0.011	ng/l	↔	,	Ω	Yes
alpha-BHC	0.011	l/gn	⊣	,	⊃	Yes
beta-BHC	0.011	ng/l	Н	ı	n	Yes
delta-BHC	0.011	l/gn	ш	•	Ω	Yes
gamma-BHC (Lindane)	0.011	l/gn	₩	,	⊃	Yes
alpha-Chlordane	0.011	l/gn	1	•	⊃	Yes
gamma-Chlordane	0.011	l/gn	1	1	⊃	Yes
Dieldrin	0.011	l/gn	1	•	⊃	Yes
4,4'-DDD	0.011	l/gn	~	•	n	Yes
4,4'-DDE	0.011	ng/l	1	•	⊃	Yes
4,4'-DDT	0.011	ug/i	П	•	⊃	Yes
Endrin	0.011	ng/l	1	ι	J	Yes
Endosulfan sulfate	0.011	l/Bn		,	⊃	Yes
Endrin aldehyde	0.011	ng/l	-	1)	Yes
Endrin ketone	0.011	ng/l	₩	1	⊃	Yes
Endosulfan-I	0.011	ng/l	₩	ı	⊃	Yes
Endosulfan-II	0.011	ng/l	₽	•	⊃	Yes
Heptachlor	0.011	l/Bn	₩	ı	n	Yes
Heptachlor epoxide	0.011	l/8n	₽	,	⊃	Yes
Methoxychlor	0.022	ng/l	7	•	⊃	Yes
Toxaphene	0.27	ng/l	-	٠	⊃	Yes

Sample ID: JC36372-6 Sample location: BMSMC, Humacao, PR

Sampling date: 26-Jan-17 Matrix: AQ - Equipment Blank Soil

Analyte Name	Result	Units Dil	Dilution Factor	Lab Flag	Validation	Reportable	
Aldrin	0.011	ug/l	Т	9.)	Yes	
alpha-BHC	0.011	ng/l	⊣	.0	n	Yes	
beta-BHC	0.011	ng/l	-	a	n	Yes	
delta-BHC	0.011	ng/l	₽	c)	Yes	
gamma-BHC (Lindane)	0.011	ng/l	H	.1	n	Yes	
alpha-Chlordane	0.011	ng/l	⊣	e;	n	Yes	
gamma-Chlordane	0.011	ng/l	₩	19	⊃	Yes	
Dieldrin	0.011	ng/l	₩	,	n	Yes	
4,4'-DDD	0.011	ng/l	-	39)	Yes	
4,4'-DDE	0.011	ng/l	1	£	n	Yes	
4,4'-DDT	0.011	ng/l	-	şa	n	Yes	
Endrin	0.011	ng/l	1	1	n	Yes	
Endosulfan sulfate	0.011	ng/l	1	1	⊃	Yes	
Endrin aldehyde	0.011	ng/l	1	e.	⊃	Yes	
Endrin ketone	0.011	l∕8n	1	1	n	Yes	
Endosulfan-I	0.011	l∕8n	T		n	Yes	
Endosulfan-II	0.011	ng/l	1	1	D	Yes	
Heptachlor	0.011	l∕8n	1	×	⊃	Yes	
Heptachlor epoxide	0.011	ng/l	1	6	n	Yes	
Methoxychlor	0.022	ng/l	1	1	⊃	Yes	
Toxaphene	0.27	l/gn	₽		D	Yes	

Sample ID: JC36372-4MS

Sample location: BMSMC, Humacao, PR Sampling date: 26-Jan-17 Matrix: Soil

	1					
Analyte Name	Result	Units Dil	Jnits Dilution Factor	Lab Flag	Lab Flag Validation	Reportable
Aldrin	46.9	ug/kg	1	1	,	Yes
alpha-BHC	17.4	ug/kg	Н	ı	•	Yes
beta-BHC	0.6	ug/kg	Н	1	,	Yes
delta-BHC	19.5	ug/kg	П	1	,	Yes
gamma-BHC (Lindane)	18.5	ug/kg	Н	ı	•	Yes
alpha-Chlordane	17.0	ug/kg	П	1	•	Yes
gamma-Chlordane	12.3	ug/kg	Н	•	-	Yes
Dieldrin	16.9	ug/kg	1	•	1	Yes
4,4'-DDD	17.2	ug/kg	1	•	•	Yes
4,4'-DDE	12.9	ug/kg	1	•	-,	Yes
4,4'-DDT	17.8	ug/kg	1	•	•	Yes
Endrin	22.6	ug/kg	1	•	,	Yes
Endosulfan sulfate	14.3	ug/kg	Н	•	,	Yes
Endrin aldehyde	25.8	ug/kg	Н	ŧ	,	Yes
Endosulfan-I	16.3	ug/kg	1	1	•	Yes
Endosulfan-il	16.4	ug/kg	1	•	•	Yes
Heptachlor	19.9	ug/kg	1	•	•	Yes
Heptachlor epoxide	18.0	ug/kg	1	•		Yes
Methoxychlor	20.3	ug/kg	1	•		Yes
Endrin ketone	19.6	ug/kg	⊣		,	Yes
Toxaphene	QN	ug/kg	П	1	•	Yes

Sample ID: JC36372-4MSD

. . . .

Sample location: BMSMC, Humacao, PR Sampling date: 26-Jan-17 Matrix: Soil

	1						
Analyte Name	Result	Units Dil	Jnits Dilution Factor	Lab Flag	Lab Flag Validation	Reportable	
Aldrin	17.0	ug/kg	1	•	•	Yes	
alpha-BHC	17.5	ug/kg	П	ı	•	Yes	
beta-BHC	16.9	ug/kg	н	1	•	Yes	
delta-BHC	19.8	ug/kg	Н	1	1	Yes	
gamma-BHC (Lindane)	18.1	ug/kg	Н	1	•	Yes	
alpha-Chlordane	20.1	ug/kg	н	•	•	Yes	
gamma-Chlordane	16.4	ug/kg	П	1	•	Yes	
Dieldrin	17.7	ug/kg	П	,	•	Yes	
4,4'-DDD	16.6	ug/kg	н	1	•	Yes	
4,4'-DDE	16.4	ug/kg	1	•		Yes	
4,4'-DDT	16.7	ug/kg	П	1	•	Yes	
Endrin	21.2	ug/kg	₩	•	1	Yes	
Endosulfan sulfate	14.8	ug/kg	Т	•	1	Yes	
Endrin aldehyde	16.3	ug/kg		1	,	Yes	
Endosulfan-l	16.7	ug/kg	1		,	Yes	
Endosulfan-II	16.4	ug/kg	Н	•	•	Yes	
Heptachlor	19.7	ug/kg	ᆏ	,	,	Yes	
Heptachlor epoxide	18.9	ug/kg	Н	1	•	Yes	
Methoxychlor	15.9	ug/kg	Н	1	•	Yes	
Endrin ketone	16.2	ug/kg	П	k	•	Yes	
Toxaphene	QN	ug/kg	H	ř	•:	Yes	

	Project/CasNumber:JC36372 Sampling Date:01/26/2017 Shipping Date:01/26/17 EPA Region No.:2
REVIEW OF PESTICIDE ORG	GANIC PACKAGE
The following guidelines for evaluating volatile organization actions. This document will assist the remake more informed decision and in better servers ample results were assessed according to USEP/IIII the following order of precedence Hazardous William Nation 1, June, 2015. SOM02.2. Pesticide Data validation actions listed on the data review workdocument, unless otherwise noted.	viewer in using professional judgment to ving the needs of the data users. The A data validation guidance documents in Vaste Support Section SOP No. HW-36A Validation. The QC criteria and data
The hardcopied (laboratory name) _Accutest reviewed and the quality control and performance data summa	data package received has beer arized. The data review for VOCs included:
Lab. Project/SDG No.:JC36372 No. of Samples:8 Trip blank No.:	X Laboratory Control SpikesX Field DuplicatesX CalibrationsX Compound IdentificationsX Compound QuantitationX Quantitation Limits
1	ound not detected ted nondetect

DATA COMPLETENESS

MISSING INFORMATION	DATE LAB. CONTACTED	DATE RECEIVED
1		
	^	
-		
<u> </u>		<u></u>
-		
		\ <u></u>
		*
		-

All criteria were met _	_X	_
Criteria were not met		
and/or see below		

HOLDING TIMES

The objective of this parameter is to ascertain the validity of the results based on the holding time of the sample from time of collection to the time of analysis.

Complete table for all samples and note the analysis and/or preservation not within criteria

SAMPLE ID	DATE	DATE	ACTION
	SAMPLED	EXTRACTED/ANALYZED	<u> </u>
Samples properly pro	eserved. All sample	es extracted and analyzed wit	thin the required criteria.
		·	

Note:

<u>Criteria</u>

Aqueous samples - seven (7) days from sample collection for extraction; 40 days from sample collection for analysis.

Non-aqueous samples – fourteen (14) days from sample collection for extraction; 40 days from sample collection for analysis.

Cooler temperature (Criteria: 4 ± 2 °C): 3.6°C - OK

<u>Actions</u>

Qualify aqueous sample results using preservation and technical holding time information as follows:

- a. If there is no evidence that the samples were properly preserved (T = 4° C \pm 2° C), and the samples were extracted or analyzed within the technical holding times, qualify detects as estimated (J) and non-detects as estimated (UJ).
- b. If there is no evidence that the samples were properly preserved ($T = 4^{\circ}C \pm 2^{\circ}C$), and the samples were extracted or analyzed outside the technical holding times, qualify detects as estimated (J) and non-detects as estimated (UJ).
- c. If the samples were properly preserved, and were extracted and analyzed within the technical holding times, no qualification of the data is necessary.
- d. If the samples were properly preserved, and were extracted or analyzed outside the technical holding times, qualify detects as estimated (J) and non-detects as estimated (UJ). Note in the Data Review Narrative that holding times were exceeded and the effect of exceeding the holding time on the resulting data.

- e. Use professional judgment to qualify samples whose temperature upon receipt at the laboratory is either below 2 degrees centigrade or above 6 degrees centigrade.
- f. If technical holding times are grossly exceeded, use professional judgment to qualify the data.

Qualify non-aqueous sample results using preservation and technical holding time information as follows:

- a. If there is no evidence that the samples were properly preserved (T = 4° C \pm 2° C), and the samples were extracted or analyzed within the technical holding time, qualify detects as estimated (J) and non-detects as estimated (UJ).
- b. If there is no evidence that the samples were properly preserved (T = 4° C \pm 2° C), and the samples were extracted or analyzed outside the technical holding time, qualify detects as estimated (J) and non-detects as estimated (UJ).
- c. If the samples were properly preserved, and were extracted and analyzed within the technical holding time, no qualification of the data is necessary.
- d. If the samples were properly preserved, and were extracted or analyzed outside the technical holding time, qualify detects as estimated (J) and non-detects as estimated (UJ). Note in the Data Review Narrative that holding times were exceeded and the effect of exceeding the holding time on the resulting data.
- e. Use professional judgment to qualify samples whose temperature upon receipt at the laboratory is either below 2 degrees centigrade or above 6 degrees centigrade.
- f. If technical holding times are grossly exceeded, use professional judgment to qualify the data.

	All criteria were met	X
Criteria	were not met see below	

GAS CHROMATOGRAPH WITH ELECTRON CAPTURE DETECTOR (GC/ECD) INSTRUMENT PERFORMANCE CHECK (SECTIONS 1 TO 5)

1. Resolution Check Mixture

Criteria

Is the resolution between two adjacent peaks in the Resolution Check Mixture C greater than or equal to 80.0% for all analytes for the primary column and greater than or equal to 50.0% for the confirmation column?

Yes? or No?

Is the resolution between two adjacent peaks in the Resolution Check Mixture (A and B) greater than or equal to 60.0%?

Yes? or No?

Note: If resolution criteria are not met, the quantitative results may not be accurate due to inadequate resolution. Qualitative identifications may also be questionable if coelution exists.

Action

- a. Qualify detects for target compounds that were not adequately resolved as tentatively identified (NJ).
- b. Qualify non-detected compounds as unusable (R).

2. Performance Evaluation Mixture (PEM) Resolution Criteria

Criteria

Is PEM analysis performed at the required frequency (at the end of each pesticide initial calibration sequence and every 12 hours)?

Yes? or No?

Action

a. If PEM is not performed at the required frequency, qualify all associated sample and blank results as unusable (R).

Criteria

Is PEM % Resolution < 90%?

Yes? or No?

Action

- a. a. Qualify detects for target compounds that were not adequately resolved as tentatively identified (NJ).
- b. Qualify non-detected compounds as unusable (R).

	All criteria were metX	
Criteria	were not met see below	

3. PEM 4,4'-DDT Breakdown

Criteria

Is the PEM 4,4'-DDT % Breakdown >20.0% and 4,4'-DDT is detected?

Yes? or No?

Action

a. Qualify detects for 4,4'-DDT; detects for 4,4'-DDD; and detects for 4,4'-DDE as estimated (J)

Criteria

Is the PEM 4,4'-DDT % Breakdown >20.0% and 4,4'-DDT is not detected

Yes? or No?

Action

- a. Qualify non-detects for 4,4'- DDT as unusable (R)
- b. Qualify detects for 4,4'-DDD as tentatively identified (NJ)
- c. Qualify detects for 4,4'-DDE as tentatively identified (NJ)

4. PEM Endrin Breakdown

Criteria

Is the PEM Endrin % Breakdown >20.0% and Endrin is detected?

Yes? or No?

Action

a. Qualify detects for Endrin; detects for Endrin aldehyde; and detects for Endrin ketone as estimated
 (J)

Criteria

Is the PEM Endrin % Breakdown >20.0% and Endrin is not detected

Yes? or No?

Action

- a. Qualify non-detects for Endrin as unusable (R)
- b. Qualify detects for Endrin aldehyde as tentatively identified (NJ)
- c. Qualify detects for Endrin ketone as tentatively identified (NJ)

All criteria were met_	_X	
Criteria were not met see belo	w	

5. Mid-point Individual Standard Mixture Resolution -

Criteria

Is the resolution between two adjacent peaks in the Resolution Check Mixture C greater than or equal to 80.0% for all analytes for the primary column and greater than or equal to 50.0% for the confirmation column?

Yes? or No?

Is the resolution between two adjacent peaks in the Resolution Check Mixture (A and B) greater than or equal to 90.0%?

Yes? or No?

Note: If resolution criteria are not met, the quantitative results may not be accurate due to inadequate resolution. Qualitative identifications may also be questionable if coelution exists.

Action

- a. Qualify detects for target compounds that were not adequately resolved as tentatively identified (NJ).
- b. Qualify non-detected compounds as unusable (R).

Criteria

Is mid-point individual standard mixture analysis performed at the required frequency (every 12 hours)?

Yes? or No?

Action

a. If the mid-point individual standard mixture analysis is not performed at the required frequency, qualify all associated sample and blank results as unusable (R).

All criteria were met	_X	
Criteria were not met		
and/or see below		_

CALIBRATION VERIFICATION

Compliance requirements for satisfactory instrument calibration are established to ensure that the instrument is capable of producing and maintaining acceptable quantitative data.

Date of initial calibration:	01/23/17	11/15/16
Dates of initial calibration verificat		11/15/16
Dates of continuing calibration:	_02/01/17:_02/02/17;_02/07/17	02/02/17
Dates of final calibration		02/02/17
Instrument ID numbers:	GC4G	GC8H
Matrix/Level:	Aqueous/low	Aqueous/low
	•	

DATE	LAB ID#	FILE	CRITERIA OUT RFs, %RSD, %D, r	COMPOUND	SAMPLES AFFECTED

Note: Initial and initial calibration verification within the guidance document performance criteria. Continuing calibration % differences meet the performance criteria in at least one of the two columns.

Final calibration verification not included in data package. No action taken.

Criteria

Are a five point calibration curve delivered with concentration levels as shown in Table 3 of SOP HW-36A, Revision 0, June, 2015?

Yes? or No?

Actions

If the standard concentrations listed in Table 3 are not used, use professional judgment to evaluate the effect on the data

Criteria

Are RT Windows calculated correctly?

Yes? or No?

Action

Recalculate the windows and use the corrected values for all evaluations.

Criteria

Are the Percent Relative Standard Deviation (%RSD) of the CFs for each of the single component target compounds less than or equal to 20.0%, except for alpha-BHC and delta-BHC?

Yes? or No?

All criteria were met_	X
Criteria were not met	
and/or see below	1.0

Are the %RSD of the CFs for alpha-BHC and delta-BHC less than or equal to 25.0%. Yes? or No?

Is the %RSD of the CFs for each of the Toxaphene peaks must be < 30% when 5-point ICAL is performed?

Yes? or No?

Is the %RSD of the CFs for the two surrogates (tetrachloro-m-xylene and decachlorobiphenyl) less than or equal to 30.0%. Yes? or No?

Action

- a. If the %RSD criteria are not met, qualify detects as estimated (J) and use professional judgment to qualify non-detected target compounds.
- b. If the %RSD criteria are within allowable limits, no qualification of the data is necessary

Continuing Calibration Checks

Criteria

Is the continuing calibration standard analyzed at the acceptable time intervals? Yes? or No?

Action

- a. If more than 14 hours has elapsed from the injection of the instrument blank that begins an analytical sequence (opening CCV) and the injection of either a PEM or mid-point concentration of the Individual Standard Mixtures (A and B) or (C), qualify all data as unusable (R).
- b. If more than 12 hours has elapsed from the injection of the instrument blank that begins an analytical sequence (opening CCV) and the injection of the last sample or blank that is part of the same analytical sequence, qualify all data as unusable (R).
- c. If more than 72 hours has elapsed from the injection of the sample with a Toxaphene detection and the Toxaphene Calibration Verification Standard (CS3), qualify all data as unusable (R).

Criteria

Is the Percent Difference (%D) within ±25.0% for the PEM sample?

Yes? or No?

Action

a. Qualify associated detects as estimated (J) and non-detects as estimated (UJ).

Criteria

For the Calibration Verification Standard (CS3); is the Percent Difference (%D) within ± 25.0%? Yes? or No?

Action

Qualify associated detects as estimated (J) and non-detects as estimated (UJ).

Criteria

Is the PEM 4,4'-DDT % Breakdown >20.0% and 4,4'-DDT is detected?

Yes? or No?

Action

- a. Qualify detects for 4,4'-DDT; detects for 4,4'-DDD; and detects for 4,4'-DDE as estimated (J)
- b. Non-detected associated compounds are not qualified

Criteria

Is the PEM 4,4'-DDT % Breakdown >20.0% and 4,4'-DDT is not detected

Yes? or No?

Action

- a. Qualify non-detects for 4,4'- DDT as unusable (R)
- b. Qualify detects for 4,4'-DDD as tentatively identified (NJ)
- c. Qualify detects for 4,4'-DDE as tentatively identified (NJ)

Criteria

Is the PEM Endrin % Breakdown >20.0% and Endrin is detected?

Yes? or No?

Action

- a. Qualify detects for Endrin; detects for Endrin aldehyde; and detects for Endrin ketone as estimated (J)
- b. Non-detected associated compounds are not qualified

Criteria

Is the PEM Endrin % Breakdown >20.0% and Endrin is not detected

Yes? or No?

Action

- a. Qualify non-detects for Endrin as unusable (R)
- b. Qualify detects for Endrin aldehyde as tentatively identified (NJ)
- c. Qualify detects for Endrin ketone as tentatively identified (NJ)

All criteria were metX
Criteria were not met
and/or see below

BLANK ANALYSIS RESULTS (Sections 1 & 2)

The assessment of the blank analysis results is to determine the existence and magnitude of contamination problems. The criteria for evaluation of blanks apply only to blanks associated with the samples, including trip, equipment, and laboratory blanks. If problems with any blanks exist, all data associated with the case must be carefully evaluated to determine whether or not there is an inherent variability in the data for the case, or if the problem is an isolated occurrence not affecting other data.

CRQL concentrat	ionN/	A		
Laboratory blanks	;			
DATE ANALYZED	LAB ID	LEVEL/ MATRIX	COMPOUND	CONCENTRATION UNITS
ug/L				it_of_0.01,_0.02,_and_0.25_
Field/Equipment	/Trip blank			
			COMPOUND	CONCENTRATION
DATE ANALYZED	LAB ID	LEVEL/ MATRIX	COMPOUND	CONCENTRATION UNITS
ANALYZED		MATRIX		
ANALYZED _No_target_anay	te_detected_in	MATRIX _the_field/equ	ipment_blanks_analyzed	UNITS d_with_this_data_package
ANALYZED _No_target_anay	te_detected_in	MATRIX _the_field/equ	ipment_blanks_analyzed	UNITS d_with_this_data_package
ANALYZED _No_target_anay	le_detected_in	MATRIX _the_field/equ	ipment_blanks_analyzed	UNITS d_with_this_data_package
ANALYZED _No_target_anayl	te_detected_in	MATRIX _the_field/equ	ipment_blanks_analyzed	UNITS d_with_this_data_package
ANALYZED _No_target_anayl	te_detected_in	MATRIX _the_field/equ	ipment_blanks_analyzed	UNITS d_with_this_data_package

All criteria were met __X__ Criteria were not met

BLANK ANALYSIS RESULTS (Section 3)

Blank Actions

Action Levels (ALs) should be based upon the highest concentration of contaminant determined in any blank. Do not qualify any blank with another blank. The ALs for samples which have been diluted should be corrected for the sample dilution factor and/or % moisture, where applicable. No positive sample results should be reported unless the concentration of the compound in the samples exceeds the ALs:

The concentration of non-target compounds in all blanks must be less than or equal to 10 μ g/L. The concentration of each target compound found in the method or field blanks must be less than its CRQL listed in the method.

Data concerning the field blanks are not evaluated as part of the CCS process. If field blanks are present, the data reviewer should evaluate this data in a similar fashion as the method blanks.

Specific actions are as follows:

Blank Actions for Pesticide Analyses

Blank Type	Blank Result	Sample Result	Action for Samples
	Detects	Not detected	No qualification required
	< CRQL	< CRQL	Report CRQL value with a U
		≥ CRQL	No qualification required
Method, Sulfur		< CRQL	Report CRQL value with a U
Cleanup, Instrument, Field, TCLP/SPLP	> CRQL	≥ CRQL and ≤ blank concentration	Report blank value for sample concentration with a U
		≥ CRQL and > blank concentration	No qualification required
	= CRQL	≤ CRQL	Report CRQL value with a U
		> CRQL	No qualification required
	Gross contamination	Detects	Report blank value for sample concentration with a U

All criteria were met __X__

Criteria	were	not	met	
and/or s	see be	elow		

CONTAMINATION SOURCE/LEVEL	COMPOUND	CONC/UNITS	AL/UNITS	SQL	AFFECTED SAMPLES
				 	- X Y X

All criteria were met __X__ Criteria were not met and/or see below _____

SURROGATE SPIKE RECOVERIES

Laboratory performance of individual samples is established by evaluation of surrogate spike recoveries. All samples are spiked with surrogate compounds prior to sample analysis. The accuracy of the analysis is measured by the surrogate percent recovery. Since the effects of the sample matrix are frequently outside the control of the laboratory and may present relatively unique problems, the validation of data is frequently subjective and demands analytical experience and professional judgment.

List the percent recoveries (%Rs) which do not meet the criteria for surrogate recovery.

Matrix:_Aqueous/Solid					_	
Lab	Lab					
Sample ID	File ID	S1 a	S1 b	\$2 a	S2 b	
JC36372-5	8G2635.D	90	89	56	48	
JC36372-6	8G2636.D	86	85	45	38	
OP233-BS1	8G2577.D	87	76	102	80	
OP233-MB1	8G2576.D	88	81	71	61	
OP233-MB1	4G77729.D	96	92	71	70	
OP233-MS	8G2631.D	83	85	106	90	
OP233-MSD	8G2632.D	85	85	102	87	
JC36372-1	4G77547.D	84	78	130	86	
JC36372-2	4G77548.D	88	81	128	80	
JC36372-3	4G77549.D	89	82	126	87	
JC36372-4	4G77550.D	84	78	118	73	
OP215-BS1	4G77546.D	91	85	122	97	
OP215-MB1	4G77545.D	88	81	125	92	
OP215-MS	4G77551.D	83	76	112	78	
OP215-MSD	4G77552.D	88	80	103	69	
Surrogate Compounds		Recovery Limits (Aqueous)			ous)	Recovery limits (Solid
S1 = Tetrachloro-m-xyl	ene	13-153	13-153%			24-136 %
S2 = Decachlorobiphenyl		10-138%				10-153 %

⁽a) Recovery from GC signal #1

Note: Surrogate recoveries were within laboratory control limits.

Actions:

⁽b) Recovery from GC signal #2

⁽c) Outside control limits due to matrix interference with the internal standard.

- a. For any surrogate recovery greater than 150%, qualify detected target compounds as biased high (J+).
- b. Do not qualify non-detected target compounds for surrogate recovery > 150 %.
- c. If both surrogate recoveries are greater than or equal to 30% and less than or equal to 150%, no qualification of the data is necessary.
- d. For any surrogate recovery greater than or equal to 10% and less than 30%, qualify detected target compounds as biased low (J-).
- e. For any surrogate recovery greater than or equal to 10% and less than 30%, qualify non-detected target compounds as approximated (UJ).
- f. If low surrogate recoveries are from sample dilution, professional judgment should be used to determine if the resulting data should be qualified. If sample dilution is not a factor:
 - i. Qualify detected target compounds as biased low (J-).
 - ii. Qualify non-detected target compounds as unusable (R).
- g. If surrogate RTs in PEMs, Individual Standard Mixtures, samples, and blanks are outside of the RT Windows, the reviewer must use professional judgment to qualify data.
- h. If surrogate RTs are within RT windows, no qualification of the data is necessary.
- i. If the two surrogates were not added to all samples, MS/MSDs, standards, LCSs, and blanks, use professional judgment in qualifying data as missing surrogate analyte may not directly apply to target analytes.

Summary Surrogate Actions for Pesticide Analyses

	Action*			
Criteria	Detected Target Compounds	Non-detected Target Compounds		
%R > 150%	J+	No qualification		
30% < %R < 150%	No qua	alification		
10% < %R < 30%	J-	UJ		
%R < 10% (sample dilution not a factor)	J-	R		
%R < 10% (sample dilution is a factor)	Use professional judgment			
RT out of RT window	Use professional judgment			
RT within RT window	No qua	alification		

* Use professional judgment in qualifying data, as surrogate recovery problems may not directly apply to target analytes.

All criteria were met _	
Criteria were not met	
and/or see below	_X

Method: SW846 8081B

MATRIX SPIKE/MATRIX SPIKE DUPLICATE (MS/MSD)

This data is generated to determine long term precision and accuracy in the analytical method for various matrices. This data alone cannot be used to evaluate the precision and accuracy of individual samples. If any % R in the MS or MSD falls outside the designated range, the reviewer should determine if there are matrix effects, i.e. LCS data are within the QC limits but MS/MSD data are outside QC limit.

1. MS/MSD Recoveries and Precision Criteria

Data for MS and MSDs will not be present unless requested by the Region.

Notify the Contract Laboratory Program Project Officer (CLP PO) if a field blank was used for the MS and MSD, unless designated as such by the Region.

NOTE: For a Matrix Spike that does not meet criteria, apply the action to only the field sample used to prepare the Matrix Spike sample. If it is clearly stated in the data validation materials that the samples were taken through incremental sampling or some other method guaranteeing the homogeneity of the sample group, then the entire sample group may be qualified.

List the %Rs, RPD of the compounds which do not meet the criteria.

Sample ID: ___JC36372-4MS/MSD ___ Matrix/Level: __Soil ____ Sample ID: ___JC36371-2MS/MSD ___ Matrix/Level: _Aqueous ___

The QC reported here applies to the following samples: JC36372-1, JC36372-2, JC36372-3, JC36372-4

Compound	JC36372-4 ug/kg Q	Spike ug/kg	MS ug/kg	MS %	Spike ug/kg	MSD ug/kg	MSD %	RPD	Limits Rec/RPD
Aldrin	ND	19	46.9	247* a	19	17.0	90	94* a	23-143/44
beta-BHC	ND	19	9.0	47	19	16.9	89	61* a	7-143/48

⁽a) Outside the QC limits.

Note: MS/MSD % recoveries and RPD within laboratory control limits except for the cases described in this document. Results for Aldrin qualified as estimated (J) in sample JC36372-4. No qualifications made based on RPD results, professional judgment.

Action

No qualification of the data is necessary on MS and MSD data alone. However, using professional judgment, the validator may use the MS and MSD results in conjunction with other QC criteria and determine the need for some qualification of the data.

A separate worksheet should be used for each MS/MSD pair.

^{*} Outside QC limits

All criteria were met _	_X_	_
Criteria were not met		
and/or see below		_

LABORATORY CONTROL SAMPLE (LCS) ANALYSIS

This data is generated to determine accuracy of the analytical method for various matrices.

1. LCS Recoveries Criteria

LCS Spike Compound	Recovery Limits (%)
gamma-BHC	50 – 120
Heptachlor epoxide	50 – 150
Dieldrin	30 – 130
4,4'-DDE	50 – 150
Endrin	50 – 120
Endosulfan sulfate	50 – 120
trans-Chlordane	30 – 130
Tetrachloro-m-xylene (surrogate)	30 – 150
Decachlorobiphenyl (surrogate)	30 – 150

COMPOUND	% R	
		QC LIMIT
y_and_RPD_within_laborat	tory_control_limits	

	y_and_RPD_within_laborat	y_and_RPD_within_laboratory_control_limits

Action

The following guidance is suggested for qualifying sample data for which the associated LCS does not meet the required criteria.

- a. If the LCS recovery exceeds the upper acceptance limit, qualify detected target compounds as estimated (J). Do not qualify non-detected target compounds.
- b. If the LCS recovery is less than the lower acceptance limit, qualify detected target compounds as estimated (J) and non-detects as unusable (R).
- c. Use professional judgment to qualify data for compounds other than those compounds that are included in the LCS.
- d. Use professional judgment to qualify non-LCS compounds. Take into account the compound class, compound recovery efficiency, analytical problems associated with each compound, and comparability in the performance of the LCS compound to the non-LCS compound.
- e. If the LCS recovery is within allowable limits, no qualification of the data is necessary.

2. Frequency Criteria:

Where LCS analyzed at the required frequency and for each matrix? <u>Yes</u> or No. If no, the data may be affected. Use professional judgment to determine the severity of the effect and qualify data accordingly. Discuss any actions below and list the samples affected.

All criteria were met	
Criteria were not met	
and/or see belowN/A	

FLORISIL CARTRIDGE PERFORMANCE CHECK

NOTE: Florisil cartridge cleanup is mandatory for all extracts.

Criteria

Is the Florisil cartridge performance check conducted at least once on each lot of cartridges used for sample cleanup or every 6 months, whichever is most frequent? Yes? or No? N/A

Criteria

Are the results for the Florisil Cartridge Performance Check solution included with the data package?

Yes? or No? N/A

Note: If % criteria are not met, examine the raw data for the presence of polar interferences and use professional judgment in qualifying the data as follows:

Action:

- a. If the Percent Recovery is greater than 120% for any of the pesticide target compounds in the Florisil Cartridge Performance Check, qualify detected compounds as estimated (J). Do not qualify non-detected target compounds.
- b. If the Percent Recovery is greater than or equal to 80% and less than or equal to 120% for all the pesticide target compounds, no qualification of the data is necessary.
- c. If the Percent Recovery is greater than or equal to 10% and less than 80% for any of the pesticide target compounds in the Florisil Cartridge Performance Check, qualify detected target compounds as estimated (J) and non-detected target compounds as approximated (UJ).
- d. If the Percent Recovery is less than 10% for any of the pesticide target compounds in the Florisil Cartridge Performance Check, qualify detected compounds as estimated (J) and qualify non-detected target compounds as unusable (R).
- e. If the Percent Recovery of 2,4,5-trichlorophenol in the Florisil Cartridge Performance Check is greater than or equal to 5%, use professional judgment to qualify detected and non-detected target compounds, considering interference on the sample chromatogram.

Note: State in the Data Review Narrative potential effects on the sample data resulting from the Florisil Cartridge Performance Check analysis not yielding acceptable results.

Note:_ No information for Florisil cartridge performance check included in data package.

There is evidence that Florisil cartridge was used for sample extraction/clean-up. No qualification of the data performed, professional judgment.

All criteria were met	
Criteria were not met	
and/or see below	

GEL PERMEATION CHROMATOGRAPHY (GPC) PERFORMANCE CHECK

NOTE: GPC cleanup is mandatory for all soil samples.

If GPC criteria are not met, examine the raw data for the presence of high molecular weight contaminants; examine subsequent sample data for unusual peaks; and use professional judgment in qualifying the data. Notify the Contract Laboratory Program Project Officer (CLP PO) if the laboratory chooses to analyze samples under unacceptable GPC criteria.

Action:

- a. If the Percent Recovery is less than 10% for the pesticide compounds and surrogates during the GPC calibration check, the non-detected target compounds may be suspect, qualify detected compounds as estimated (J).
- b. If the Percent Recovery is less than 10% for the pesticide compounds and surrogates during the GPC calibration check, qualify all non-detected target compounds as unusable (R).
- c. If the Percent Recovery is greater than or equal to 10% and is less than 80% for any of the pesticide target compounds in the GPC calibration, qualify detected target compounds as estimated (J) and non-detected target compounds as approximated (UJ).
- d. If the Percent Recovery is greater than or equal to 80% and less than or equal to 120% for all the pesticide target compounds, no qualification of the data is necessary.
- e. If high recoveries (i.e., greater than 120%) were obtained for the pesticides and surrogates during the GPC calibration check, qualify detected compounds as estimated (J). Do not qualify non-detected target compounds.

Note: State in the Data Review Narrative potential effects on the sample data resulting from the GPC cleanup analyses not yielding acceptable results.

Note: No information for performance of GPC cleanup included in data package. No qualification of the data performed, professional judgment.

All criteria were met _	_X
Criteria were not met	
and/or see below	

TARGET COMPOUND IDENTIFICATION

Criteria:

- 1. Is Retention Times (RTs) of both of the surrogates and reported target compounds in each sample within the calculated RT Windows on both columns?

 Yes? or No?
- 2. Is the Tetrachloro-m-xylene (TCX) RT ±0.05 minutes of the Mean RT (RT) determined from the initial calibration and Decachlorobiphenyl (DCB) within ±0.10 minutes of the RT determined from the initial calibration?

 Yes? or No?
- 3. Is the Percent Difference (%D) for the detected mean concentrations of a pesticide target compound between the two Gas Chromatograph (GC) columns within the inclusive range of \pm 25.0 %?

 Yes? or No?
- 4. When no analytes are identified in a sample; are the chromatograms from the analyses of the sample extract and the low-point standard of the initial calibration associated with those analyses on the same scaling factor?

 Yes? or No?
- 5. Does the chromatograms display the Single Component Pesticides (SCPs) detected in the sample and the largest peak of any multi-component analyte detected in the sample at less than full scale.

 Yes? or No?
- 6. If an extract is diluted; does the chromatogram display SCPs peaks between 10-100% of full scale, and multi-component analytes between 25-100% of full scale? Yes? or No? N/A
- 7. For any sample; does the baseline of the chromatogram return to below 50% of full scale before the elution time of alpha-BHC, and also return to below 25% of full scale after the elution time of alpha-BHC and before the elution time of DCB?

 Yes? or No?
- 8. If a chromatogram is replotted electronically to meet these requirements; is the scaling factor used displayed on the chromatogram, and both the initial chromatogram and the replotted chromatogram submitted in the data package.

 Yes? or No?

Action:

- a. If the qualitative criteria for both columns were not met, all target compounds that are reported as detected should be considered non-detected.
- b. Use professional judgment to assign an appropriate quantitation limit using the following guidance:
 - i. If the detected target compound peak was sufficiently outside the pesticide RT Window, the reported values may be a false positive and should be replaced with the sample Contract Required Quantitation Limits (CRQL) value.

- ii. If the detected target compound peak poses an interference with potential detection of another target peak, the reported value should be considered and qualified as unusable (R).
- c. If the data reviewer identifies a peak in both GC column analyses that falls within the appropriate RT Windows, but was reported as a non-detect, the compound may be a false negative. Use professional judgment to decide if the compound should be included.

Note: State in the Data Review Narrative all conclusions made regarding target compound identification.

- d. If the Toxaphene peak RT windows determined from the calibration overlap with SCPs or chromatographic interferences, use professional judgment to qualify the data.
- e. If target compounds were detected on both GC columns, and the Percent Difference between the two results is greater than 25.0%, consider the potential for coelution and use professional judgment to decide whether a much larger concentration obtained on one column versus the other indicates the presence of an interfering compound. If an interfering compound is indicated, use professional judgment to determine how best to report, and if necessary, qualify the data according to these guidelines.
- f. If Toxaphene exhibits a marginal pattern-matching quality, use professional judgment to establish whether the differences are due to environmental "weathering" (i.e., degradation of the earlier eluting peaks relative to the later eluting peaks). If the presence of Toxaphene is strongly suggested, report results as presumptively present (N).

GAS CHROMATOGRAPH/MASS SPECTROMETER (GC/MS) CONFIRMATION

NOTE: This confirmation is not usually provided by the laboratory. In cases where it is provided, use professional judgment to determine if data qualified with "C" can be salvaged if it was previously qualified as unusable (R).

Action:

- a. If the quantitative criteria for both columns were met (\geq 5.0 ng/ μ L for SCPs and \geq 125 ng/ μ L for Toxaphene), determine whether GC/MS confirmation was performed. If it was performed, qualify the data using the following guidance:
 - i. If GC/MS confirmation was not required because the quantitative criteria for both columns was not met, but it was still performed, use professional judgment when evaluating the data to decide whether the detect should be qualified with "C".
 - ii. If GC/MS confirmation was performed, but unsuccessful for a target compound detected by GC/ECD analysis, qualify those detects as "X".

All criteria were met	_X
Criteria were not met	
and/or see below	

COMPOUND QUANTITATION AND REPORTED CONTRACT REQUIRED QUANTITATION LIMITS (CRQLS)

The sample quantitation evaluation is to verify laboratory quantitation results. In the space below, please show a minimum of one sample calculation:

Action:

- a. If sample quantitation is different from the reported value, qualify result as unusable (R).
- b. When a sample is analyzed at more than one dilution, the lowest CRQLs are used unless a QC exceedance dictates the use of the higher CRQLs from the diluted sample.
- c. Replace concentrations that exceed the calibration range in the original analysis by crossing out the "E" and its corresponding value on the original reporting form and substituting the data from the diluted sample.
- d. Results between the MDL and CRQL should be qualified as estimated (J).
- e. Results less than the MDL should be reported at the CRQL and qualified (U). MDLs themselves are not reported.
- f. For non-aqueous samples, if the percent moisture is less than 70.0%, no qualification of the data is necessary. If the percent moisture is greater than or equal to 70.0% and less than 90.0%, qualify detects as estimated (J) and non-detects as approximated (UJ). If the percent moisture is greater than or equal to 90.0%, qualify detects as estimated (J) and non-detects as unusable (R) (see Table).

Percent Moisture Actions for Pesticide Analysis for Non-Aqueous Samples

Criteria	Action				
	Detected Associated Non-detected Associated Compounds				
% Moisture < 70.0	No qualification				
70.0 < % Moisture < 90.0	J	UJ			
% Moisture > 90.0	J	R			

List samples which have < 50 % solids

 	_	
. = -		

Note: If any discrepancies are found, the Region's designated representative may contact the laboratory to obtain additional information that could resolve any differences. If a discrepancy remains unresolved, the reviewer must use professional judgment to decide which value is the most accurate. Under these circumstances, the reviewer may determine that qualification of data is warranted. Note in the Data Review Narrative a description of the reasons for data qualification and the qualification that is applied to the data.

Dilution performed

SAMPLE ID	DILUTION FACTOR	REASON FOR DILUTION
	-	
		-
		2

All criteria were met	_X
Criteria were not met	
and/or see below	

Matrix:

Soil

FIELD DUPLICATE PRECISION

Sample IDs:

NOTE: In the absence of QAPP guidance for validating data from field duplicates, the following action will be taken.

Field duplicates samples may be taken and analyzed as an indication of overall precision. These analyses measure both field and lab precision; therefore, the results may have more variability than laboratory duplicates which only laboratory performance. It is also expected that soil duplicate results will have a greater variance than water matrices due to difficulties associated with collecting identical field duplicate samples. Identify which samples within the data package are field duplicates. Estimate the relative percent difference (RPD) between the values for each compound. If large RPDs (> 50%) is observed, confirm identification of samples and note difference in the executive summary.

oumpio iso	00	.000.12 2,00000.12	<u> </u>	11161			
COMPOUND	SQL	SAMPLE	DUPLICATE	RPD	ACTION		
	ug/L	CONC.	CONC.				
Field duplicate analyzed with this data package. RPD within the required criteria of < 50 %.							

JC36372-2/JC36372-3

Actions:

- a. Qualify as estimated positive results (J) and nondetects (UJ) for the compound that exceeded the above criteria. For organics, only the sample and duplicate will be qualified.
- b. If an RPD cannot be calculated because one or both of the sample results is not detected, the following actions apply:
 - If one sample result is not detected and the other is greater than 5x the SQL qualify (J/UJ).
 - ii. If one sample value is not detected and the other is greater than 5x the SQL and the SQLs for the sample and duplicate are significantly different, use professional judgment to determine if qualification is appropriate.
 - iii. If one sample value is not detected and the other is less than 5x, use professional judgment to determine if qualification is appropriate.
 - iv. If both sample and duplicate results are not detected, no action is needed.

OVERALL ASSESSMENT OF DATA Action:

- 1. Use professional judgment to determine if there is any need to qualify data which were not qualified based on the Quality Control (QC) criteria previously discussed.
- 2. Write a brief narrative to give the user an indication of the analytical limitations of the data.

Note: The Contract Laboratory Program Project Officer (CLP PO) must be informed if any inconsistency of the data with the Sample Delivery Group (SDG) Narrative. If sufficient information on the intended use and required quality of the data is available, the reviewer should include their assessment of the usability of the data within the given context. This may be used as part of a formal Data Quality Assessment (DQA).

Overall assessment of the data: Results are valid; the data can be used for decision

making purposes.

EXECUTIVE NARRATIVE

SDG No:

JC36372

Laboratory:

Accutest, New Jersey

Analysis:

SW846-8015C

Number of Samples:

7

Location:

BMSMC, Humacao, PR

SUMMARY:

Seven (7) samples were analyzed for the trace metals following method SW846-6010C and SW846-7470A/7471B for Hg. The sample results were assessed according to USEPA data validation guidance documents in the following order of precedence: Hazardous Waste Support Section SOP NO. HW-3b Revision 0 (July 2015) ISM02 ICP-MS Data Validation; USEPA Contract Laboratory program National Functional Guidelines for Inorganic data Review (OSWER 9240.1-45, EPA 540-R-04-004, October 2004- Final). Validation of Metal for the Contract Laboratory Program (CLP) (SOP HW-2, Revision 13. Based on ILM05.3 (August 2009). Quality control validation criteria were derived from "Test Methods for Evaluating Solid Waste, Physical/Chemical Methods SW-846 (Final Update IV, 1998)". The QC criteria and data validation actions listed on the data review worksheets are from the primary guidance document, unless otherwise noted.

Results are valid and can be used for decision making purposes.

Critical issues:

None

Major:

None

Minor:

None

Critical findings:

None

Major findings:

None

Minor findings:

- 1. Target analytes not detected in the field/equipment blank except for the cases described in this document:
- Ba, Ca, Mn and Na detected in the field blank at a concentration below the reporting limit. The laboratory qualified the results with a B qualifier. No further
- Ba, Ca, Mn, Hg and Na detected in the equipment blank at a concentration below the reporting limit.

The laboratory qualified the results with a B qualifier. No further qualification performed.

2. MS/MSD % recovery and RPD within laboratory control limits except for the cases described in the Data Review Worksheet.

No action taken on samples with high analyte concentration compared to the amount spiked.

No action taken on samples with either the MS or MSD % recovery within laboratory and guidance document control limits, professional judgment

No qualification made based on RPD outside laboratory control limits, professional judgment.

Results for Sb qualified as estimated (J or UJ) in all samples.

No MS/MSD sample analyzed for the aqueous matrix, except for Hg.

3. Positive results with concentration between the MDL and the reporting limit (RL) are qualified as estimated (J or UJ).

COMMENTS:

Results are valid and can be used for decision making purposes.

Reviewers Name:

Rafael Infante

Chemist License 1888

Rafael Defaut

Signature:

Date:

February 18, 2017

SAMPLE METAL DATA SAMPLE SUMMARY

Sample ID: JC36372-1

Sample location: BMSMC, Humacao, PR

Sampling date: 26-Jan-17

Matrix: Soil

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Aluminum	21800	ug/kg	1.0	•	- 9	Yes
Antimony	0.39	ug/kg	1.0	В	JB ✓	Yes
Arsenic	2.0	ug/kg	1.0	В	JB	Yes
Barium	113	ug/kg	1.0	-	- 1	Yes
Beryllium	0.21	ug/kg	1.0	В	JB V	Yes
Cadmium	0.14	ug/kg	1.0	В	JB 🗸	Yes
Calcium	6510	ug/kg	1.0	-	-	Yes
Chromium	11.9	ug/kg	1.0	-	-	Yes
Cobalt	15.2	ug/kg	1.0	-	-	Yes
Copper	39.2	ug/kg	1.0	-	-	Yes
Iron	35000	ug/kg	2.0	-	-	Yes
Lead	7.7	ug/kg	1.0	-	-	Yes
Magnesium	7630	ug/kg	1.0	-	-	Yes
Manganese	873	ug/kg	1.0	-	- //	Yes
Mercury	0.031	ug/kg	1.0	В	JB V	Yes
Nickel	5.5	ug/kg	1.0	-	1	Yes
Potassium	1580	ug/kg	1.0	-	-	Yes
Selenium	0.51	ug/kg	1.0	U	U	Yes
Silver	0.22	ug/kg	2.0	U	U	/ Yes
Sodium	353	ug/kg	1.0	В	JB 🗸 🗸	Yes
Thallium	0.44	ug/kg	1.0	U	U	Yes
Vanadium	94.3	ug/kg	1.0	-	-	Yes
Zinc	139	ug/kg	1.0	-	-	Yes

Sample ID: JC36372-2

Sample location: BMSMC, Humacao, PR

Sampling date: 26-Jan-17

Matrix: Soil

	•					
Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Aluminum	18000	ug/kg	1.0	-	-	Yes
Antimony	0.63	ug/kg	1.0	В	JB √	✓ Yes
Arsenic	1.6	ug/kg	1.0	В	JB	Yes
Barium	98.9	ug/kg	1.0	-	I	Yes
Beryllium	0.16	ug/kg	1.0	В	JB /	Yes
Cadmium	0.16	ug/kg	1.0	В	JB 🗸	Yes
Calcium	10000	ug/kg	1.0	-	-	Yes
Chromium	20.4	ug/kg	1.0	-	-	Yes
Cobalt	13.3	ug/kg	1.0	-	-	Yes
Copper	38.3	ug/kg	1.0	-	-	Yes
Iron	29200	ug/kg	1.0	-	-	Yes
Lead	9.2	ug/kg	1.0	•	-	Yes
Magnesium	7230	ug/kg	1.0	-	-	Yes
Manganese	797	ug/kg	1.0	-	-	Yes
Mercury	0.050	ug/kg	1.0	-	-	Yes
Nickel	13.5	ug/kg	1.0	-	-	Yes
Potassium	2210	ug/kg	1.0	-	-	Yes
Selenium	0.58	ug/kg	1.0	U	U	Yes
Silver	0.12	ug/kg	1.0	U	U	Yes
Sodium	144	ug/kg	1.0	В	JB (/	Yes
Thallium	0.50	ug/kg	1.0	U	U	Yes
Vanadium	81.8	ug/kg	1.0	-	-	Yes
Zinc	129	ug/kg	1.0	-	-	Yes

Sample ID: JC36372-4

Sample location: BMSMC, Humacao, PR

Sampling date: 26-Jan-17

Matrix: Soil

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Aluminum	17000	ug/kg	1.0	-	-	Yes
Antimony	0.70	ug/kg	1.0	В	JB V	Yes
Arsenic	1.8	ug/kg	1.0	В	JB 🗸 🗸	Yes
Barium	97.5	ug/kg	1.0		- ,	Yes
Beryllium	0.16	ug/kg	1.0	В	JB 🗸	Yes
Cadmium	0.17	ug/kg	1.0	В	JB 🗸 🗸	Yes
Calcium	5110	ug/kg	1.0	5.		Yes
Chromium	12.6	ug/kg	1.0	-	-	Yes
Cobalt	10.8	ug/kg	1.0	*	•	Yes
Copper	35.3	ug/kg	1.0	2:		Yes
Iron	27200	ug/kg	1.0	-	-	Yes
Lead	29	ug/kg	1.0	2	2	Yes
Magnesium	4820	ug/kg	1.0	*		Yes
Manganese	621	ug/kg	1.0		•	Yes
Mercury	0.10	ug/kg	1.0	-	-	Yes
Nickel	6.2	ug/kg	1.0	-	-	Yes
Potassium	1180	ug/kg	1.0		-	Yes
Selenium	0.51	ug/kg	1.0	U	U	Yes
Silver	0.11	ug/kg	1.0	U	U //	Yes
Sodium	183	ug/kg	1.0	В	JB V	Yes
Thallium	0.44	ug/kg	1.0	U	U	Yes
Vanadium	81.3	ug/kg	1.0	2	-	Yes
Zinc	85.1	ug/kg	1.0	27	-	Yes

Sample ID: JC36372-4MS

Sample location: BMSMC, Humacao, PR

Sampling date: 26-Jan-17

Matrix: Soil

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Aluminum	31800	ug/kg	1.0	-	-	Yes
Antimony	103	ug/kg	1.0	-	-	Yes
Arsenic	247	ug/kg	1.0	•	-	Yes
Barium	437	ug/kg	1.0	-	-	Yes
Beryllium	257	ug/kg	1.0	-	-	Yes
Cadmium	252	ug/kg	1.0	-	-	Yes
Calcium	11600	ug/kg	1.0	-	-	Yes
Chromium	272	ug/kg	1.0	-	-	Yes
Cobalt	270	ug/kg	1.0	-	-	Yes
Copper	305	ug/kg	1.0	-	-	Yes
Iron	39000	ug/kg	1.0	-	-	Yes
Lead	287	ug/kg	1.0	-	-	Yes
Magnesium	9100	ug/kg	1.0	-	-	Yes
Manganese	1450	ug/kg	1.0	•	-	Yes
Mercury	0.373	ug/kg	1.0	-	-	Yes
Nickel	265	ug/kg	1.0	-	-	Yes
Potassium	4870	ug/kg	1.0	•	-	Yes
Selenium	247	ug/kg	1.0	-	-	Yes
Silver	30.1	ug/kg	1.0	-	-	Yes
Sodium	3550	ug/kg	1.0	-	-	Yes
Thallium	250	ug/kg	1.0	-	-	Yes
Vanadium	346	ug/kg	1.0	-	-	Yes
Zinc	339	ug/kg	1.0	-	-	Yes

Sample ID: JC36372-4MSD

Sample location: BMSMC, Humacao, PR

Sampling date: 26-Jan-17

Matrix: Soil

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Aluminum	25900	ug/kg	1.0	-	-	Yes
Antimony	103	ug/kg	1.0	-	-	Yes
Arsenic	228	ug/kg	1.0	-	-	Yes
Barium	338	ug/kg	1.0	+	-	Yes
Beryllium	237	ug/kg	1.0	70	(5)	Yes
Cadmium	234	ug/kg	1.0	-	-	Yes
Calcium	8440	ug/kg	1.0	- 	-	Yes
Chromium	245	ug/kg	1.0	-	-	Yes
Cobalt	245	ug/kg	1.0	7.0	(5.0	Yes
Copper	274	ug/kg	1.0	21	-	Yes
Iron	29800	ug/kg	1.0	5	-	Yes
Lead	256	ug/kg	1.0	2	-	Yes
Magnesium	7450	ug/kg	1.0	*		Yes
Manganese	848	ug/kg	1.0	0		Yes
Mercury	0.378	ug/kg	1.0	7.	-	Yes
Nickel	242	ug/kg	1.0	0.	2	Yes
Potassium	4280	ug/kg	1.0	-	(*)	Yes
Selenium	229	ug/kg	1.0	-	2.1	Yes
Silver	28.6	ug/kg	1.0	*	-	Yes
Sodium	3230	ug/kg	1.0	-	-	Yes
Thallium	234	ug/kg	1.0	~	-	Yes
Vanadium	309	ug/kg	1.0		-	Yes
Zinc	305	ug/kg	1.0	*	90	Yes

Sample ID: JC36372-5

Sample location: BMSMC, Humacao, PR

Sampling date: 26-Jan-17

Matrix: AQ - Field Blank Soil

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Aluminum	21	ug/l	1.0	U	U	Yes
Antimony	3.3	ug/l	1.0	U	U	Yes
Arsenic	2.2	ug/l	1.0	U	U	Yes
Barium	1.6	ug/l	1.0	В	JB 🗸	Yes
Beryllium	0.25	ug/l	1.0	U	U	Yes
Cadmium	0.4	ug/l	1.0	U	U	Yes
Calcium	188	ug/l	1.0	В	JB 💆	Yes
Chromium	0.81	ug/l	1.0	U	U	Yes
Cobalt	0.69	ug/l	1.0	U	U	Yes
Copper	2.4	ug/l	1.0	U	U	Yes
Iron	12	ug/l	1.0	U	U	Yes
Lead	2.3	ug/l	1.0	U	U	Yes
Magnesium	85	ug/l	1.0	U	U	Yes
Manganese	5.1	ug/l	1.0	В	JB	Yes
Mercury	0.047	ug/l	1.0	U	U	Yes
Nickel	0.76	ug/l	1.0	U	U	Yes
Potassium	120	ug/l	1.0	U	U	Yes
Selenium	4.1	ug/l	1.0	U	U	Yes
Silver	0.88	ug/l	1.0	U	U	Yes
Sodium	154	ug/l	1.0	В	JB 🗸	Yes
Thallium	1.9	ug/l	1.0	U	U	Yes
Vanadium	0.66	ug/l	1.0	U	U	Yes
Zinc	1.3	ug/l	1.0	U	U	Yes

Sample ID: JC36372-6

Sample location: BMSMC, Humacao, PR

Sampling date: 26-Jan-17

Matrix: AQ - Equipment Blank

Analyte Name	Result	Units	Dilution Factor	Lab Flag	Validation	Reportable
Aluminum	21	ug/l	1.0	U	U	Yes
Antimony	3.3	ug/l	1.0	U	U	Yes
Arsenic	2.2	ug/l	1.0	U	U	Yes
Barium	1.9	ug/l	1.0	В	JB 🗸	Yes
Beryllium	0.25	ug/l	1.0	U	U	Yes
Cadmium	0.4	ug/l	1.0	U	U	Yes
Calcium	86.9	ug/l	1.0	В	JB	Yes
Chromium	0.81	ug/l	1.0	U	U	Yes
Cobalt	0.69	ug/l	1.0	U	U	Yes
Copper	2.4	ug/l	1.0	U	U	Yes
Iron	12	ug/l	1.0	U	U	Yes
Lead	2.3	ug/l	1.0	U	U	Yes
Magnesium	85	ug/l	1.0	U	U	Yes
Manganese	2.1	ug/l	1.0	В	JB	Yes
Mercury	0.067	ug/l	1.0	В	JB /	Yes
Nickel	1.00	ug/l	1.0	В	JB Y	Yes
Potassium	120	ug/l	1.0	U	U	Yes
Selenium	4.1	ug/l	1.0	U	U	Yes
Silver	0.88	ug/l	1.0	U	U	Yes
Sodium	65.7	ug/l	1.0	В	JB √	Yes
Thallium	1.9	ug/l	1.0	U	U	Yes
Vanadium	0.66	ug/l	1.0	U	U	Yes
Zinc	1.3	ug/l	1.0	U	U	Yes

Type of validation	Full:X	Project Number:JC36372
Type of validation	Limited:	Date:01/26/17
	EPA Region:2_	Date:01/26/17
F	REVIEW OF INORGANIO	C ANALYSIS DATA PACKAGE
sulfide, and/or cyanidassist the reviewer in serving the needs of the validation guidance do Section SOP NO. HW-Laboratory program NA5, EPA 540-R-04-06 Program (CLP) (SOF validation criteria were Methods SW-846 (Fininformation (if available)	e were created to deline using professional judgrine data users. The sampocuments in the following above the deline of the following and the following and the following are seen as a few of the following and the following are seen as a few of the following and the following are seen as a few of the following are seen as a few of the few of	als analyses (6010C/6020/7000A series method) eate required validation actions. This document will ment to make more informed decision and in better ole results were assessed according to USEPA data gorder of precedence: Hazardous Waste Support 15) ISM02 ICP-MS Data Validation; USEPA Contract elines for Inorganic data Review (OSWER 9240.1-17). Validation of Metal for the Contract Laboratory Based on ILM05.3 (August 2009). Quality controphods for Evaluating Solid Waste, Physical/Chemical The project QAPP is reviewed for project specific data validation actions listed on the data review sument, unless otherwise noted.
		est data package received has been been been been been been been bee
No. of Samples: Field blank No.:	JC36372-6	Sample matrix:Soil
X Data deliver:X Holding TimeX CalibrationsX BlanksX ICP Interfere	ables es	X Laboratory DuplicatesX_ Field DuplicatesX_ Laboratory Control SamplesX_ ICP Serial Dilution ResultsX_ Detection Limits ResultsX_ Sample Quantitation
Overall Comments: _M	Metal_Analysis_(SW846-	-6010C);_Hg_(SW846-7470A/7471B)
Definition of Qualifiers		
J- Estimated res U- Compound no R- Rejected data UJ- Estimated nor E- Laboratory qui Reviewer:	ot detected	
Date:02/18/2017		

			All criteria were metX Criteria were not met and/or see below
1.	DATA DELIVERABLI	ES	
	A. Data Packag	e :	
MISSI	NG INFORMATION	DATE LAB. CONTACTED	DATE RECEIVED
		3772.22	
_			N 12 / / / / / / / / / / / / / / / / / /
	B. Other Discre	pancies:	
-			
		2.14.66. 22.81	
4 370		TO 3 TO 10 10 10 10 10 10 10 10 10 10 10 10 10	
	-		<u> </u>

All criteria were metX	
Criteria were not met	
and/or see below	

HOLDING TIMES

The objective of this parameter is to ascertain the validity of the results based on the holding time of the sample from time of collection to the time of preparation, and subsequently from the time of preparation to the time of analysis.

Complete table for all samples and circle the analysis date for samples not within criteria

SAMPLE ID	DATE SAMPLED	CYANIDE DATE ANALYSIS	Hg DATE ANALYSIS	OTHERS DATE ANALYSIS	рH	SULFIDE	ACTION
SAMPLES [DIGESTED AN	D ANALYZE	D WITHIN T	HE METHO	D REC	OMMEND	ED HOLDING

<u>Criteria</u>

Metals – 180 days from time of collection.
Mercury – 28 days from time of collection.
Hexavalent Chromium (solids)- 30/7 from day of collection; 48 hrs aqueous samples
Cyanide – 14 days from time of collection
Sulfide - 14 days from time of collection
pH measurements of aqueous samples upon receipt at the laboratory (criteria pH ≤ 2 for metals
pH ≥ 12 for cyanide)

Actions: Qualify positive results/nondetects as follows:

If holding times are exceeded, estimate positive results (J) and rejects nondetects (R).
If pH > 2 for metals or pH < 12 for cyanide, positive results (J) and nondetects (UJ).
Cooler Temperature (Criteria: 4°C + 2°C):3.6°C
If cooler temperature is > 10°C, flag non-detects as (UJ) and detects as (J).

All criteria were metN/A
Criteria were not met
and/or see below

ICP-MS TUNE ANALYSIS

Is the ICP-MS tuned prior to calibration?

Yes or No?

Does the % RSD exceeds 5% for any isotope in the tuning solution?

Yes or No?

Action:

NOTES: For ICP-MS tunes that do not meet the technical criteria, apply the action to all samples reported from the analytical run.

- 1. If the ICP-MS instrument was not tuned prior to calibration, the sample data should be qualified as unusable (R).
- 2. If the tuning solution was not analyzed or scanned at least 5x consecutively or the tuning solution does not contain the required analytes spanning the analytical range, the reviewer should use professional judgment to determine if the associated sample data should be qualified. The reviewer may need to obtain additional information from the laboratory. The situation should be recorded in the Data Review Narrative and noted for Contract Laboratory Program Project Officer (CLP PO) action.
- 3. If the resolution of the mass calibration is not within 0.1 u for any isotope in the tuning solution, qualify all analyte results that are ≥ Method Detection Limit (MDL) associated with that isotope as estimated (J), and all non-detects associated with that isotope as estimated (UJ). The situation should be recorded in the Data Review Narrative and noted for CLP PO action.
- 4. If the %RSD exceeds 5% for any isotope in the tuning solution, qualify all sample results that are ≥ MDL associated with that tune as estimated (J), and all non-detects associated with that tune as estimated (UJ). The situation should be recorded in the Data Review Narrative and noted for CLP PO action.

Table 2. ICP-MS Tune Actions for ICP-MS Analysis

ICP-MS Tune Results	Action for Samples		
Tune not performed	Qualify all results as unusable (R)		
Tune not performed properly	Use professional judgment		
Resolution of mass calibration not within 0.1u	Qualify results that are ≥ MDL as estimated (J)		
	Qualify non-detects as estimated (UJ)		
% RSD > 5%	Qualify results that are ≥ MDL as estimated (J)		
	Qualify non-detects as estimated (UJ)		

Note:

All criteria were metX
Criteria were not met
and/or see below

INSTRUMENT CALIBRATION (SECTION 1)

Compliance requirements for satisfactory instrument calibration are established to ensure that the instrument is capable of producing and maintaining acceptable quantitative data. Minimum of 2 calibration points for ICP-AES and ICP-MS; 5 points for Hg; and 4 points for cyanide. One initial calibration standard at the CRQL level for cyanide and Hg. If no, write in the non-compliance section of the data review narrative.

List the analytes which did not meet the percent recovery (%R) criteria for Initial or Continuing Calibration Verification standards (ICV or CCV).

Acceptance Criteria	ICV %R	CCV %R
Metals by 6010C/6020	100 + 10%	100 + 10%
Mercury/Metals by 7000s	100 + 10%	100 + 20%
Cyanide	100 + 15%	100 + 15%
Sulfide	100 + 15%	100 + 15%

DATE	ICV/CCV#	ANALYTE	%R	ACTION	SAMPLES AFFECTED
INITI	AL AND CONTI	NUING CALIBRA	ATION M	LEET METHOD SPEC	IFIC CRITERIA
		· · · · · · · · · · · · · · · · · · ·			

ACTIONS: If any analyte does not meet the %R criteria, follow the actions stated below. Qualify five samples on either side of the ICV/CCV out of control limit.

Estimate positive results (J) if: Metals by 6010C/6020 Mercury/Metals by 7000s Cyanide Sulfide	ICV 111 – 125% 111 – 125% 116 – 130% 116 – 130%	CCV 111 – 125% 111 – 135% 116 – 130% 116 – 130%
Estimate positive results and nondetects (U/Us Metals by 6010C/6020 Mercury/Metals by 7000s Cyanide Sulfide	J) if: 75 – 89% 75 – 89% 70 – 84% 70 – 84%	75 – 89% 65 – 79% 70 – 84% 70 – 84%
Reject positive results and nondetects (R) if: Metals by 6010C/6020 Mercury/Metals by 7000s Cyanide Sulfide	<75%, >125% <75%, >125% <70%, >130% <70%, >130%	<75%, >125% <65%, >135% <70%, >130% <70%, >130%

All criteria were met __X___ Criteria were not met and/or see below

- III. INSTRUMENT CALIBRATIONS (SECTIONS 2 & 3)
- Analytical Sequence

Did the laboratory use the proper number of standards for calibration as described in the method?

Yes or No

B. Were calibrations performed at the beginning of each analysis?

Yes or No

Were calibration verification standards analyzed at the beginning of sample analysis and the proper frequency according to the method?

Yes_or No

Where the AA correlation coefficients (r) for the calibration curves
 ≥ 0.995? If r < 0.995, estimate positive results and nondetects (J/UJ).
 It is not necessary to qualify results if the laboratory used order regression.

Yes or No

Data quality may be affected if any of the above answer are "no". Use professional judgment to determine the severity of the effect and qualify the data accordingly. Discuss any actions below and list the sample affected.

Other Check Standards

Laboratories may analyze an additional check standard after establishing the calibration curve. This standard may contain low level concentrations of target analytes and be analyzed and evaluated by the laboratory similar to a CLP "CRLD" standard (CRI for ICP, CRA for AA, and/or mid-range standard for CN and Sulfide). A $100 \pm 20\%$ recovery acceptance limit should be used by the validator to evaluate the standard.

ACTIONS: If any analyte does not meet the %R criteria, follow the action needed below. Qualify 50% of either side of the CRI/CRA out of control limits.

% R		%R < 50%	%R 79%	=	50-	%R 150%	=	121-	%R 150%	>	Affecte	ed Ra	nge
Qualify Positi	ve/N	ondetects Res	ults										
Metals 6010C/6020	by	R/R	J/UJ			J/A			R/A		<2x CI	RI cor	nc.
Hg/metals 7000s	by	R/R	J/UJ			J/A			R/A		<1.5x conc.		CRI
Cyanide		R/R	J/UJ			J/A			R/A		<1.5x conc.	mid	std.
Sulfide		R/R	J/UJ		i.	J/A			R/A		<1.5x conc.	mid	std.

CRI is not required for Al, Ba, Ca, Fe, Mg, Na, and K.

NOTE: CRLD standard within laboratory and method specific criteria.

All criteria were met	N/A
Criteria were	not met
and/or see below	

Table 4. Calibration Actions for ICP-MS Analysis

Calibration Result	Action for Samples
Calibration not performed	Qualify all results as unusable (R)
Calibration incomplete	Use professional judgment
·	Qualify results that are ≥ MDL as estimated
	(J)
	Qualify non-detects as estimated (UJ)
Not at least one calibration standard at or	Qualify results that are ≥ MDL but < 2x the
below the CRQL for each analyte	CRQL as estimated (J)
	Qualify non-detects as estimated (UJ)
Correlation coefficient < 0.995; %D outside	Qualify results that are ≥ MDL as estimated
±30%; y-intercept ≥ CRQL	(J)
	Qualify non-detects as estimated (UJ)
Correlation coefficient < 0.990	Qualify results that are ≥ MDL as estimated
	(J)
	Qualify non-detects as unusable (R)
ICV/CCV %R < 75%	Qualify results that are ≥ MDL as unusable
	(R)
	Qualify all non-detects as unusable (R)
ICV/CCV %R 75-89%	Qualify results that are ≥ MDL as estimated
	low (J-)
	Qualify non-detects as estimated (UJ)
ICV/CCV %R 111-125%	Qualify results that are ≥ MDL as estimated
	high (J+)
ICV/CCV %R > 125%	Qualify results that are ≥ MDL as estimated
	high (J+)
ICV/CCV %R > 160%	Qualify results that are ≥ MDL as unusable
	(R)

All criteria were met	_
Criteria were not me	ŧ
and/or see belowX	_

IV. BLANK ANALYSIS RESULTS (Sections 1 & 2)

The assessment of the blank analysis results is to determine the existence and magnitude of contamination problems. The criteria for evaluation of blanks apply only to blanks associated with the samples, including equipment, field, and laboratory blanks. If problems with any blanks exist, all data associated with the case must be carefully evaluated to determine whether or not there is an inherent variability in the data for the case, or if the problem is an isolated occurrence not affecting other data.

List the contamination in Sections 1 & 2 below. A separate worksheet page should be used for soil and water blanks.

Laboratory blan	ks		Matrix:A	Aqueous	
DATE ANALYZED	ICB/CCB#	PREP BLK	ANALYTE	CONCENTRATION UNITS	
No_analyte_d		hod_blanks_	_above_reporting_limits	queous	
	_	rielel D	_	CONCENTRATION	
DATE ANALYZED	EQUIPMENT BLANK	I/FIELD	ANALYTE	UNITS	
Field/equipme	ent_blank_analy	zed_as_par	t_of_this_data_package.		

Note: Terget analytes not detected in the field/equipment blank except for the cases described in this document:

- Ba, Ca, Mn and Na detected in the field blank at a concentration below the reporting limit. The laboratory qualified the results with a B qualifier. No further
- Ba, Ca, Mn, Hg and Na detected in the equipment blank at a concentration below the reporting limit.

The laboratory qualified the results with a B qualifier. No further qualification performed.

Table. Field/Rinsate/Trip Blank Actions for ICP-MS Analysis

Blank Result	Sample Result	Action for Samples	
> CRQL	≥ MDL but ≤ CRQL	Report CRQL value with a "U	
	> CRQL but < Blank Result	Report at level of Blank Result with a "U"	
	> Blank Result but < 10x the Blank Result	Use professional judgment to qualify results as estimated (J)	

Al	Il criteria were metX Criteria were not met and/or see below
IV. BLANK ANALYSIS RESULTS (Section 3)	
Frequency requirements	
Was the preparation blank analyzed for each matrix, at the frequency of the method? If no, estimate positive results < 10x IDL for which preparation blank was If more than 20 samples/batch, qualification begins at the 21 st sample.	Yes or No as not analyzed.
B. Was an ICB analyzed?	Yes or No
C. Was a CCB analyzed at the frequency stated in the method?	Yes or No
Data quality may be affected if any of the above answer is "no". Use determine the severity of the effect and qualify the data accordingly. Data the samples affected.	
NOTE FOR SOIL SAMPLES Compare raw sample value with blank results in ug/L unit, or Convert blanks analyzed during a soil case to mg/Kg in order to compresults.	pare them with the sample
Conc. In ug/L x [Volume diluted to (mL)]/[Weight digested] x 1L/1000m 1mg/1000 □g = concentration in wet weight (mg/Kg)	nL x 1000g/1Kg x
Concentration, dry weight (mg/Kg) = (Wet weight concentration)/(% Sol	ids) x 100
BLANK ANALYSIS RESULTS (Sections 4,5)	
Laboratory blanks (PB, ICB/CCB) must first be used to qualify field and samples. Any contamination remaining in the field or equipment blank will be use samples.	

			All criteria were met Criteria were no and/or see below _	t met
4. Initial/C	Continuing Calibr	ration Blanks (ICB/CCB	s) Actions	
Are all ICB/CCI	Bs less than the	SQL?	Yes or No	
		ther side of the ICB/CC he ICB/CCB value.	B out of control limits.	
ICB/CCB#	ANALYTE	CONC/UNITS	SAMPLES AFFECTED	
Are the PB less	than the SQL?		Yes or No	
If yes, reject all	results (R) < 10	x the PB value.		
РВ	ANALYTE	CONC/UNITS	SAMPLES AFFECTED	
· · · · · · · · · · · · · · · · · · ·				
BLANK ANALY	SIS RESULTS	(Section 6)		
6. Field/E	quipment Blank	(FB/EB) Actions		
Are the	FB/EB less tha	n the SQL?	Yes or No	
If no, was the F	B/EB value alre	ady rejected due to oth	er QC criteria? Yes or No	
lf no, reject (R) the FB/EB valu		≤5x the FB/EB value.	Reject soil data with raw digest results	s < 5x
РВ	ANALYTE	CONC/UNITS	SAMPLES AFFECTED	

All criteria were met __N/A__ Criteria were not met and/or see below ____

Table 5. Calibration/Preparation Blank Actions for ICP-MS Analysis - Summary

Blank Type	Blank Result	Sample Result	Action for Samples			
ICB/CCB	≥ MDL but ≤ CRQL	Non-detect	No action			
≥ MDL but ≤ CRQL		Report CRQL value with a "U"				
> CRQL	5.09	Use professional judgment				
ICB/CCB	> CRQL	≥ MDL but ≤ CRQL	Report CRQL value with a "U"			
> CRQL but < Blank Res	sult	Report at level of Blank	Result with a "U"			
> Blank Result		Use professional judgme	ent			
ICB/CCB	≤ (-MDL) but ≥ (-CRQL)	≥ MDL, or non-detect	Use professional judgment			
ICB/CCB	< (-CRQL)	< 10x the CRQL	Qualify results that are ≥ CRQL as estimated low (J-)			
			Qualify non-detects as estimated (UJ)			
Preparation Blank	> CRQL	≥ MDL but ≤ CRQL	Report CRQL value with a "U"			
> CRQL but < 10x the B	lank Result	Qualify results as estima	ated high (J+)			
≥ 10x the Blank Result		No action				
Preparation Blank	≥ MDL but ≤ CRQL	Non-detect	No action			
≥ MDL but ≤ CRQL		Report CRQL value with	a "U"			
> CRQL		Use professional judgment				
Preparation Blank	< (-CRQL)	< 10x the CRQL	Qualify results that are ≥ CRQL as estimated low (J-)			
			Qualify non-detects as estimated (UJ)			

>					Cr	were metX iteria were not met r see below
INDUCTIVELY CO	OUPLED PLAS	SMA (ICP) IN	TERFERE	NCE CHEC	K SAMPLE	
The assessment interelement and b				mple (ICS)	is to verify	y the laboratory's
1. Recovery	Criteria					
List any elements : %).	in the ICS AB	and ICS A so	olutions wh	nich did not	meet the %F	R criteria (80 – 120
DATE EL	EMENT	%R ACT	ION	SAMPLES	S AFFECTE)
_Interference_che	ck_sample_wi	thin_method	_performa	nce_criteria		
ACTIONS: If an element does	not meet the	%R criteria,	follow the a	actions state	ed below	
% R	%R < 50%	%R = 5	50- %R 150%	= 121-	%R 150%	>
Qualify Positive/No	ondetects Res				,	
Metals by 6010C/6020	R/R	J/UJ	J/A		R/A	
2. Frequency	requirements	;				
Were interference (beginning of the a			quency sta	ted in the m		es or No
If no, <u>ACTIONS:</u> Estima	te positive res	ults (J) all sa	mples for v	which Al, Ca	a, Fe, Mg > I	CS value.
The data may be a qualify the data ac						ty of the effect and ected.
	227 (137		- 272			72 - 72
	- 1					21.51

Αll	criteria	were	met	1	\/A_	
		Crit	eria	were	not	met
	а	nd/or	see	belov	V	

Table 6. Interference Check Actions for ICP-MS Analysis - Summary

Interference Check Sample Results	Action for Samples
ICS not analyzed	Qualify detects and non-detects as unusable (R)
ICS not analyzed in proper sequence	Use professional judgment.
ICS %R>150%	Use professional judgment
ICS %R > 120% (or greater than true value + 2x the CRQL)	Qualify results that are ≥ MDL as estimated high (J+)
ICS %R 80-12-%	No qualification
ICS %R 50-79% (or less than true value – 2x the CRQL)	Qualify results that are ≥ MDL as estimated low (J-)
	Qualify non-detects as estimated (UJ)
ICSAB %R < 50%	Qualify detects as estimated low (J-) and non- detects as unusable (R)
Potential false positives in field samples with interferents	Qualify results that are ≥ MDL as estimated high (J+)
Potential false negatives in field samples with interferents	Qualify results that are ≥ MDL but < 10x the (negative value) as estimated low (J-) Qualify non-detects as estimated (UJ)

All criteria were met _	
Criteria were	not met
and/or see below	X

VI. MATRIX SPIKE (MS)

Sample # _JC36372-4MS/-4MSD_____ Matrix: __Soil__ Units: __ug/kg__ Sample # _JC36275-2AMS/-2AMSD_(Hg)_ Matrix: __Aqueous__ Units: __ug/l___

This data is generated to determine long term precision and accuracy in the analytical method for various matrices. Note that for Region 2, MS not required for: Ca, Mg, K, and Na for aqueous matrix.

Al, Ca, Fe, Mg, K, Na, for soil matrix

MS Recovery Criteria. List the percent recoveries for analytes which did not meet the %R criteria (75 – 125%); (85 – 115 % FOR Cr (VI)).

ANALYTE	SPIKE SAMPLE	SAMPLE	SPIKE	% R	ACTION
	RESULT (SSR)	RESULT (SR)	ADDED		
MS/MSD	recoveries and RPD	within laboratory	control lim	its excep	ot for the cases described in
		this do	cument.		
Al	31800	17000	2880	513.3	No action
Sb	103	0.70	231	44.3	Qualify results (J or UJ)
Ba	437	97.5	231	147.2	No action^
Ca	11600	5110	2880	225.1	No action^
Fe	39000	27200	2880	409.2	No action
Mg	9100	4820	2880	148.4	No action^
Mn	1450	621	231	359.4	No action^
K	4870	1180	2880	228.0	No action^
Al	25900	17000	2880	317.9	No action
Sb	103	0.70	231	45.7	Qualify results (J or UJ)

Note: No action taken, sample concentration high compared to the amount spiked.

^- no action taken, professional judgment. Either the MS/MSD % recovery within laboratory and guidance document control limits.

No qualification made based on RPD outside laboratory control limits, professional judgment.

Results for Sb qualified as estimated (J or UJ) in all samples.

No MS/MSD sample analyzed for the aqueous matrix, except for Hg.

ACTIONS: Matrix spike actions apply to all samples of the same matrix. The qualification will also be applied to the results of all samples within a given area of the site, if deemed appropriate.

If the sample results \geq 4x the spike concentration, no action is taken. If any analyte does not meet the %R criteria, follow the actions stated below.

Table 9. Spike Sample Actions for ICP-MS Analysis

Spike Sample Results	Action for Samples
Matrix Spike %R < 30% Post-digestion spike %R < 75%	Qualify affected results that are ≥ MDL as estimated low (J-) and affected non-detects as unusable (R)
Matrix Spike %R < 30% Post-digestion spike %R ≥ 75%	Qualify affected results that are ≥ MDL as estimated (J) and affected non-detects as estimated (UJ)
Matrix Spike %R 30-74% Post-digestion Spike %R < 75%	Qualify affected results that are ≥ MDL as estimated low (J-) and affected non-detects as estimated (UJ)
Matrix Spike %R 30-74% Post-digestion spike %R ≥ 75%	Qualify affected results that are ≥ MDL as estimated (J) and affected non-detects as estimated (UJ)
Matrix Spike %R > 125% Post-digestion spike %R > 125%	Qualify affected results that are ≥ MDL as estimated high (J+)
Matrix Spike %R > 125% Post-digestion spike %R ≤ 125%	Qualify affected results that are ≥ MDL as estimated (J)

Spike Sample Results	Action for Samples			
Matrix Spike %R < 30% No post-digestion spike performed	Qualify affected results that are ≥ MDL as estimated low (J-) and affected non-detects as unusable (R)			
Matrix Spike %R 30-74% No post-digestion spike performed	Qualify affected results that are ≥ MDL as estimated low (J-) and non-detects as estimated (UJ)			
Matrix Spike %R > 125% No post-digestion spike performed	Qualify affected results that are ≥ MDL as estimated high (J+) Non-detects are not qualified			

2. Frequency Criteria

A. Was a matrix spike prepared at the frequency stated in the method (1/20)? Yes or No

If no, estimate positive results (J) for which analyte was not spiked. If more than 20 samples/batch, qualification begins at the 21st sample.

B. Was a field blank used as spiked sample?

Yes or No

If yes, estimate positive results (J) < 4x spike level added for the analyte.

A separate worksheet page should be used for each matrix spike

	,	All criteria were metN/A Criteria were not met and/or see below
VII. FIELD DUPLICATES		
Sample #:	Matrix:	Units:_ug/L

Field duplicate samples may be taken and analyzed as an indication of overall precision. Field duplicate analyses measure both field and lab precision; therefore, the results may have more variability than laboratory duplicates which measure only laboratory performance. It is also expected that soil duplicate results will have a greater variance than water matrices due to difficulties associated with collecting identical field duplicate samples.

List the concentrations and RPDs in the field duplicate pair. RPD criteria: \pm 20% for aqueous; \pm 35% for soil. For soil duplicates, if the % solids for the sample and its duplicate differ by more than 1%, report concentrations in ug/L and calculate RPD or difference for each analyte.

ANALYTE	SQL ug/L	SQL ug/Kg	SAMPLE RESULTS	DUPLICATE RESULTS	RPD	ACTION
Al					20.5	No action
Sb						
As						
Ba					25.5	No action
Ве						
Cd						
Ca					31.5	No action
Cr						
Co						
Cu						
Fe					26.7	No action
Pb						
Mg						
Mn					52.4	No action
Hg						
Ni						
K						
Se						
Ag						
Na						
TI						
V						
Zn						
Cyanide						
Cr(VI)						

Note: No field/laboratory duplicates analyzed with data set. MS/MSD % recoveries RPD used to assess precision. RPD within laboratory and generally acceptable control limits except for the cases described in this document. No qualification made based on RPD results

Field duplicate actions should be applied to only the sample and its duplicate.

All criteria were metN/A	١
Criteria were not	met
and/or see below _	

Actions: Indicates which criterion was used to evaluate precision by circling either the RPD or SQL for each element. If both sample and duplicate are nondetects, the RPD is not calculated (NC), no action is needed.

Table 8. Duplicate Sample Actions for ICP-MS Analysis

Duplicate Sample Results	Action for Samples
Aqueous: Both original sample and duplicate sample > 5x the CRQL and 20% < RPD < 100%	Qualify those results that are ≥ CRQL as estimated (J)
Aqueous: Both original sample and duplicate sample > 5x the CRQL and RPD ≥ 100%	Qualify those results that are ≥ CRQL as unusable (R)
Soil/Sediment: Both original sample and duplicate sample > 5x the CRQL and 35% < RPD < 120%	Qualify those results that are ≥ CRQL as estimated (J)
Soil/Sediment: Both original sample and duplicate sample > 5x the CRQL and RPD ≥ 120%	Qualify those results that are ≥ CRQL as unusable (R)
Original sample or duplicate sample ≤ 5x the CRQL (including non-detects) and absolute difference between sample and duplicate > CRQL	Qualify those results that are ≥ MDL as estimated (J) and non-detects as estimated (UJ)

A separate worksheet page should be used for each laboratory duplicate analysis

TI
V
Zn
Cr(VI)
Sulfide
Cyanide

						ia were metN/A Criteria were not met and/or see below
VIII. LABO	RATORY	DUPLIC	ATES (Section	า 1)		
Laboratory run duplicates samples to verify laboratory consistency and precision. They are a measure of laboratory performance. It is also expected that soil duplicate results will have a greater variance than water matrices due to difficulties associated with collecting identical field duplicate samples.						
1. Difference C	Criteria					
for soil). For s	oil duplic	cates, if th	ne % solids fo	eting the RPD or the sample a e RPD or differe	nd its duplic	9% for aqueous; <u>+</u> 35% ate differ by more than analyte.
Sample #				Matrix:		Units:
ANALYTE	SQL ug/L	SQL mg/Kg	SAMPLE RESULTS	DUPLICATE RESULTS	RPD	ACTION
Al		1113113				
Sb						
As						
Ва			_			
Be						
Cd						
Ca						
Cr						
Co			<u> </u>			
Cu						
Fe						
Pb						
Mg						
Mn	-					
Hg	-					
Ni	+		<u> </u>			
K	-					
Se	-	 	<u> </u>			
Ag	+	-			_	

Note: No laboratory duplicate analyzed with this data package

Laboratory duplicates actions should be applied to all other samples of the same matrix type. This qualification will also be applied to the results of all samples within a given area of the site, if deemed appropriate.

All criteria were met __X__ Criteria were not met and/or see below ____

<u>Actions:</u> Indicates which criterion was used to evaluate precision by circling either the RPD or SQL for each element. If both sample and duplicate are non-detects, the RPD is not calculated (NC), no action is needed.

Table 8. Field Duplicate Sample Actions for ICP-MS Analysis

Sample Type	Field Duplicate Result	Action for Samples
Aqueous	Sample and its field duplicate ≥ 5x the CRQL and RPD > 20%	Qualify sample and its duplicate as estimated (J)
	Sample and/or its field duplicate < 5x the CRQL and absolute difference > the CRQL	Qualify results > the MDL as estimated (J) Qualify non-detects as estimated (UJ)
Soil/Sediment	Sample and its field duplicate ≥ 5x the CRQL and RPD > 50%	Qualify sample and its duplicate as estimated (J)
	Sample and/or its field duplicate < 5x the CRQL and absolute difference > 2x the CRQL	Qualify results > the MDL as estimated (J)
		Qualify non-detects as estimated (UJ)

2. Frequency Criteria

A. Was a laboratory duplicate prepared at the frequency stated in the method (1/20)? Yes or No

If no, estimate positive results (J) for the analyte which duplicate was not performed. If more than 20 samples/batch, qualification begins at the 21st sample.

B. Was a field blank used for laboratory duplicate analysis? Yes or No

If yes, estimate positive results (J) for the analyte if field blank was used for duplicate analysis.

All criteria were metX
Criteria were not met
and/or see below

IX. LABORATORY CONTROL SAMPLE (LCS/LCSD)

The assessment of the LCSs is to determine both intralaboratory contamination and matrix specific precision and accuracy. Note that for Region 2, LCS is not required for aqueous Hg and Cyanide.

LCS Recoveries Criteria

A. Aqueous LCS/Solid LCS

List any LCS recoveries not within %R criteria (80 – 120%) and the samples affected.

DATE	ELEMENT	% R	ACTION	SAMPLES AFFECTED
 Recoveries	_within_laboratory_co	ntrol_limits_		
			<u>-</u>	
		7000		
74.5	W			NUS 25

ACTIONS: If analyte does not meet the %R criteria, follow the actions stated below:

Table 7. LCS Actions for ICP-MS Analysis

LCS Result	Action for Samples
%R 40-69%	Qualify results that are ≥ MDL as estimated low (J-) Qualify non-detects as estimated (UJ)
%R > 130%	Qualify results that are ≥ MDL as estimated high (J+)
%R 70-130%	No qualification
%R < 40%	Qualify results that are ≥ MDL as estimated low (J-) Qualify non-detects as unusable (R)
%R > 150%	Qualify detects as unusable (R); non- detects no qualification

All criteria were metX
Criteria were not met
and/or see below

2. Frequency Criteria

A. Was a laboratory control sample prepared at the frequency stated in the method (1/20)? **Yes** or No

If no, estimate positive results (J) for the analyte if LCS was not performed.

If more than 20 samples/batch, qualification begins at the 21st sample.

All criteria were metX
Criteria were not met
and/or see below

X. ICP SERIAL DILUTION ANALYSIS (Section 1)

The assessment of the ICP serial dilution analysis is to determine the precision of the laboratory through a 5x dilution.

1. Percent Difference (%D) Criteria:

	analys				ed for each the undiluted					
\$	Serial	dilutions	were	not	performed	for	the	following	target	analytes:
		ilutions we entrations	•		out analytical (re dilution.	esul	ts did i	not agree wi	thin 10%	difference
List the %	6Ds for	analytes v	vhich did	not me	eet the %D cr	iteria	(10%	/100%)		
Sample #	# .IC3	6372-4			Matrix	: S	oil		Units:	ua/ka

ANALYTE	IDL	50x IDL	SAMPLE RESULTS	SERIAL DILUTION	%D	ACTION
Al				·		
Sb	0.32	16	6.3	0.0	100	No action
As	0.24	12	16.3	10.3	36.8	No action
Ва						
Ве						
Çd	0.055	2.75	1.5	0.0	100	No action
Ca						
Сг						
Со						
Cu						
Fe						
Pb						
Mg						
Mn						
Hg					<u> </u>	
Ni					1	
K						
Se						
Ag						
Na						
TI	0.44	22	1.9	0.0	100	No action
V						
Zn						

Note: Serial dilution within method performance criteria.

All criteria were metX
Criteria were not met
and/or see below

ACTIONS: Actions apply to all samples of the same matrix. The qualification will also be applied to the results of all samples within a given area of the site, if deemed appropriate. Qualify only samples with raw results > 50x MDL.

Flag results with an (E) for elements exhibiting %D > 10%. Estimate (J) positive results > 50x MDL for elements that exhibited %D > 10 but < 100.

Reject (R) positive results > 50x MDL for elements which exhibited %D ≥ 100 %.

SERIAL DILUTION ANALYSIS (Section 2)

was not performed.

2. Frequency Criteria

A. Was a serial dilution analysis prepared as required by the method? Yes or No

If no, estimate positive results ≥ 50x MDL (J) for the analyte which serial dilution analysis

B. Was a field blank used for serial dilution analysis?

Yes or No

If yes, estimate positive results \geq 50x MDL (J) for the analyte if field blank was used for serial dilution analysis.

Table 10. Serial Dilution Actions for ICP-MS Analysis

Serial Dilution Result	Action for Samples
Aqueous: Sample concentration > 50x MDL and 10% < %D < 100%	Qualify affected results whose raw data are > MDL as estimated (J)
Aqueous: Sample concentration > 50x MDL and %D ≥ 100%	Qualify affected results whose raw data are > MDL as unusable (R)
Soil/Sediment: Sample concentration > 50x MDL and 15% < %D < 120%	Qualify affected results whose raw data are > MDL as estimated (J)
Soil/Sediment: Sample concentration > 50x MDL and %D ≥ 120%	Qualify affected results whose raw data are > MDL as unusable (R)
Interferences present	Use professional judgment

A separate worksheet page should be used for each serial dilution analysis.

	erni on	Criteria were not met and/or see below
XI.	ICP-MS INTERNAL STANDARDS	
	Are internal standard added to the sample?	Yes_or No?
	Are the proper number of internal standard added to the sample?	Yes or No?
	Is the % Relative Intensities for all internal standards in a sample response in the calibration blank?	e is within 60-125% of the Yes or No?
	Note:_ICP-OES_internal_standards_used;_relative_intensities_v _document_performance_criteria	vithin_the_guidance_

All criteria were met

NI/A

Action:

NOTE: Apply the action to the affected analytes for each sample that does not meet the internal standard criteria.

- 1. If no internal standards were analyzed with the run, the sample data should be qualified as unusable (R). Record this in the Data Review Narrative and note for CLP Project Officer (CLP PO) action.
- 2. If less than five of the required internal standards were analyzed with the run, or a target analyte(s) is (are) not associated to an internal standard, the sample data, or analyte data not associated to an internal standard should be qualified as unusable (R). Record this in the Data Review Narrative and note for CLP PO action.
- 3. If the % Relative Intensities for all internal standards in a sample is within 60-125% of the response in the calibration blank, the sample data should not be qualified.
- 4. If the %RI for an internal standard in a sample is not within the 60-125% limit, qualify the data for those analytes associated with the internal standard(s) outside the limit as follows:
 - a. If the sample was reanalyzed at a two-fold dilution with internal standard %RI within the limits, report the result of the diluted analysis without qualification. If the %RI of the diluted analysis was not within the 60-125% limit, report the results of the original undiluted analyses and qualify the data for all analytes that are ≥ Method Detection Limit (MDL) in the sample associated with the internal standard as estimated (UJ).
 - b. If the sample was not reanalyzed at a two-fold dilution, the reviewer should use professional judgment to determine the reliability of the data. The reviewer may determine that the results are estimated (J) or unusable (R).

Table 11. Internal Standard Actions for ICP-MS Analysis

Internal Standard Results	Action for Samples
No internal standards	Qualify all results as unusable (R)
< 5 of the required internal standards	Qualify all results as unusable (R)
Target analyte not associated with internal standard	Qualify all analyte results not associated with an internal standard as unusable (R)
% RI < 60% or > 125%, original sample	Do not qualify the data
reanalyzed at 2-fold dilution, and % RI of diluted	
sample analysis is between 60% and 125%	
% RI < 60% or > 125%, original sample	Qualify analytes associated with the failed
reanalyzed at 2-fold dilution, and % RI of diluted	internal standard that are ≥ MDL as estimated
sample analysis is outside the 60% to 125%	(J) and qualify associated non-detects as
limit	estimated (UJ)
Original sample not reanalyzed at 2-fold dilution	Use professional judgment
	Qualify sample results as estimated (J) or unusable ®

XII. DETECTION LIMITS RESULTS

The detection limit assessment is to verify that samples results are within instrument calibration range or linear range (ICP).

Instrument Detection Limits (IDL). Note IDL is not required for Cyanide.

- A. IDL/MDL (or lowest quantitation limit used) results were present and found to be all levels that meet the project objectives? Yes or No
- B. IDL/MDL (or lowest quantitation limit used) were not met for the following elements:
- 2. Reporting Requirements
- A. Were sample results on Form I (or equivalent) reported down to the IDL/MDL or lowest quantitation limit used for all analytes? Yes or No
- B. Were sample weights, volumes, and dilutions taken into account when reporting results (positive and nondetects)? Yes or No

If no, the reported results may be inaccurate. Request the laboratory resubmit the corrected data.

- Sediment Sample Percent Solids (% solids):
- A. Were the % solids for any sediment samples < 50% but ≥ 10%? Yes or No If yes, estimate positive results and nondetects (J/UJ) if the % solids is 10-50%. List the affected samples:______
- B. Were the % solids for any sediment samples < 10%? Yes or No If yes, reject all results (R) if the % solid is < 10%. List the affected samples: N/A
- XI. TOTAL/DISSOLVED OR INORGANIC/TOTAL ANALYTES
- A. Were any analyses performed for dissolved as well as total analytes on the same sample(s)?

 Yes or **No**
- B. Were any analyses performed for inorganic as well as total analytes on the same sample(s)?

 Yes or **No**

If yes, compare the differences between dissolved (or inorganic) and total analyte concentrations. Compute each difference as a percent of the total analyte only when both of the following conditions are fulfilled:

- (1) The dissolved (or inorganic) concentration is greater than total concentration, and
- (2) greater than or equal to 5xMDL.

1.5

				All criteria were metN/A Criteria were not met and/or see below
C.	Is any dissolved than 20%?	d (or inorganic) concentrati	on greater than Yes or <u>No</u>	its total concentration by more
D.	Is any dissolved than 50%?	d (or inorganic) concentrati	on greater than Yes or <u>No</u>	its total concentration by more
	percent differen	ce is greater than 20%, ated. If the difference is mo		dissolved/inorganic and total ject (R) both the values.
XII.	SAMPLE QUAN	NTITATION		
The sar	mple quantitation	evaluation is to verify labo	ratory quantitati	on results.
	Sample results arameters.	fall within the linear range	for ICP and wi	thin the calibration range for all
dilution	If samples resu performed?	Its were beyond the linear	range/calibratio	n range of the instrument, were
List the	affected sample	s/elements/dilution:		
In the s	pace below, plea	ase show a minimum of one	e sample calcula	ation per method:
ICP/ICF	P-MS	Computer printout		
<u>Hg/Met</u>	als by AA	Computer printout		
<u>Hexava</u>	lent Chromium			
Cyanide	<u>e</u>			
<u>Others</u>				
		llowing equation may be neoncentrations (mg/Kg):	ecessary to con-	vert raw data values reported in
Conc. ii		diluted to, mL x 1L 1000 m		ng = concentration 000 mg in wet weight

In addition the sample results are converted to dry weight by using the percent solid calculations:

Wet weight concentration x 100 = final concentration, dry weight (mg/Kg) % solids

27

OVERALL ASSESSMENT

Action:

- 1. Use professional judgment to determine if there is any need to qualify data which were not qualified based on the QC criteria previously discussed.
- 2. Write a brief Data Review Narrative to give the user an indication of the analytical limitations of the data. Note any discrepancies between the data and the Sample Delivery Group (SDG) Narrative for Contract Laboratory Program Project Officer (CLP PO) action. If sufficient information on the intended use and required quality of the data is available, the reviewer should include an assessment of the data usability within the given context.
- 3. If any discrepancies are found, the laboratory may be contacted by the Region's designated representative to obtain additional information for resolution. If a discrepancy remains unresolved, the reviewer may determine that qualification of the data is warranted.

_